PANPAC®

Panpac Medical Corporation
6F-2, No. 202, Sec. 3, Ta-Tong Road,
Shi-Chih District, New Taipei City, 22103
Taiwan
TEL: +886-2-86472242 FAX: +886-2-86472770

To be submitted to the Competent Health Authorities in Russian Federation

I herewith confirm that the content of **USER'S MANUAL** for medical device
"Urogynecological pessary" is valid and true.

Panpac Medical Corporation, Taiwan, 6F-2, No. 202, Sec.3, Ta-Tong Road, Shi-Chih Dist., New Taipei City, Taiwan, is the original manufacturer of the device mentioned in the document.

Signature: Evelyn Hsu
Wenlin Hsu 2021/6/15

(Sales Manager: Wenlin Hsu (Evelyn):

(Notary marks):

On this 18 day of June 2021 before me, the undersigned notary public, personally appeared

HSU, WENLIN, proved to me through satisfactory evidence of identification,
which was passport to be person whose name is signed on the preceding or attached
document, and acknowledged to me that he/ she signed it voluntarily for its stated purpose.

Notary Public (name) Yufen Wu

(Notary stamp):

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案號：北院民認育字 0381	號日期：2021-06-18
Case No.:	Date:
本文件之簽名或蓋章，在臺灣臺北地方法院所屬民間 公證人吳育芬事務所認證。公證人：吳育芬	
Attested at the Wu, Yu-Fen Notary Public Office of Taiwan Taipei District Court, R.O.C., that the signature(s)/seal(s) in this document is/are authentic.	
Notary Public	<u>Yufen Wu</u>

The Wu, Yu-Fen Notary Public Office,
Taiwan Taipei District Court, Republic Of China
Add: 8F - 7, No.179, Fuxing N. Rd.,
Songsshan Dist., Taipei City, Taiwan(R.O.C.)

1243-2
Panpac Medical Corporation
6F-2, No. 202, Sec. 3, Ta-Tong Road,
Shi-Chih District, New Taipei City, 22103
Taiwan
TEL: +886-2-86472242 FAX: +886-2-86472770

To be submitted to the Competent Health Authorities in Russian Federation

I herewith confirm that the content of **USER'S MANUAL** for medical device
"Urogynecological pessary" is valid and true.

Panpac Medical Corporation, Taiwan, 6F-2, No. 202, Sec.3, Ta-Tong Road, Shi-Chih Dist., New Taipei City, Taiwan, is the original manufacturer of the device mentioned in the document.

Evelyn Hsu
Signature: Wenlin Hsu 2021/6/15

(Sales Manager: Wenlin Hsu (Evelyn):

(Notary marks):

On this 18 day of June 2021 before me, the undersigned notary public, personally appeared

HSU, WENLIN, proved to me through satisfactory evidence of identification,

which was passport to be person whose name is signed on the preceding or attached document, and acknowledged to me that he/ she signed it voluntarily for its stated purpose.

Notary Public (name) *yu fen*

(Notary stamp):

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案號：北院民認育字 0381	號日期：2021-06-18
Case No.:	Date:
本文件之簽名或蓋章，在臺灣臺北地方法院所屬民間公證人吳育芬事務所認證。公證人：吳育芬	
Attested at the Wu, Yu-Fen Notary Public Office of Taiwan Taipei District Court, R.O.C., that the signature(s)/seal(s) in this document is/are authentic.	
Notary Public	<i>yu fen</i> Wu, Yu-Fen

The Wu, Yu-Fen Notary Public Office,
Taiwan Taipei District Court, Republic Of China
Add: 8F - 7, No.179, Fuxing N. Rd.,
Songshan Dist., Taipei City, Taiwan(R.O.C.)

PANPAC®

Panpac Medical Corporation
6F-2, No. 202, Sec. 3, Ta-Tong Road,
Shi-Chih District, New Taipei City, 22103
Taiwan
TEL: +886-2-86472242 FAX: +886-2-86472770

APPROVED BY

Evelyn Hsu

Sales Manager

Evelyn Hsu

Wenlin Hsu 2021/6/15

Sales Manager

Wenlin Hsu (Evelyn)

USER'S MANUAL

of medical device

"Urogynecological pessary"

manufactured by Panpac Medical Corporation, Taiwan, 6F-2, No. 202, Sec.3, Ta-Tong Road,
Shi-Chih Dist., New Taipei City, Taiwan, 22103.

V. 03

June 2021

Table of contents

1. Name of medical device.....	3
2.1. Manufacturer:.....	5
2.2. Authorized Representative in the Russian Federation.....	5
3. Appointment.....	5
3.1. Scope.....	5
3.2. Contraindications.....	6
3.3. Possible side effects.....	6
3.4. Precautionary measures.....	6
4. Performance characteristics, therapeutic indications.....	7
5. Application risks.....	9
6. Specifications.....	10
6.1. Silicone hardness.....	26
7. Accessories Description.....	26
8. Drug content.....	26
9. Mode of application.....	26
9.1. Commissioning.....	26
9.2. Cleaning, disinfection, sterilization before first use.....	27
9.3. Selection of product size.....	28
9.4. Pessary installation.....	28
9.5. Female Patient Instruction.....	32
10. Instructions for reprocessing and sterilization.....	32
11. Pessary extraction.....	33
12. Duration of pessaries use.....	33
13. Circumstances requiring medical advice.....	33
14. Terms of Use.....	33
15. Information contained in the leaflet.....	33
16. Material Information.....	34
17. Environmental requirements.....	34
19. Transportation.....	35
20. Storage conditions, shelf life, service life.....	35
21. Packaging and labeling of a medical device.....	35
21.1. Individual packing.....	35
21.2. Individual packaging labeling.....	36
21.3. Individual packaging labeling project for the Russian Federation.....	37
21.4. Labeling layout for the Russian Federation.....	38
22. Disposal.....	39
23. Maintenance and Repair Requirements.....	39
24. Warranty.....	39
25. Complaints.....	40

1. Name of medical device

**"Urogynecological pessary",
design variants**

1. Pessary Maryland (Marland) without support consisting of:
 - 1.1. Pessary, size: 2, 3, 4, 5, 6, 7, 8,
 - 1.2. Pessary extractor PRM2 (if required),
 - 1.3. Insert instruction.
2. Pessary Maryland (Marland) with support from:
 - 2.1. Pessary, size: 2, 3, 4, 5, 6, 7, 8,
 - 2.2. Pessary extractor PRM2 (if required),
 - 2.3. Insert instruction.
3. Shaatz pessary with drainage consisting of:
 - 3.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8,
 - 3.2. Pessary extractor PRM2 (if required),
 - 3.3. Insert instruction.
4. Pessary Cup (Cup) without support consisting of:
 - 4.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8,
 - 4.2. Pessary extractor PRM2 (if required),
 - 4.3. Insert instruction.
5. Pessary Cup (Cup) with support from:
 - 5.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7,
 - 5.2. Pessary extractor PRM2 (if required),
 - 5.3. Insert instruction.
6. Pessary Gellhorn with short stem and drainage consisting of:
 - 6.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8,
 - 6.2. Pessary extractor PRM2 (if required),
 - 6.3. Insert instruction.
7. Pessary Gellhorn with a long stem and drainage consisting of:
 - 7.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8,
 - 7.2. Pessary extractor PRM2 (if required),
 - 7.3. Insert instruction.
8. Pessary Gehrung consisting of:
 - 8.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8,
 - 8.2. Pessary extractor PRM2 (if required),
 - 8.3. Insert instruction.
9. Pessary Ring (Ring) without support consisting of:
 - 9.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,
 - 9.2. Pessary extractor PRM2 (if required),
 - 9.3. Insert instruction.
10. Pessary Ring (Ring) with support consisting of:
 - 10.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,
 - 10.2. Pessary extractor PRM2 (if required),
 - 10.3. Insert instruction.
11. Ring w / Knob pessary without support consisting of:
 - 11.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7,
 - 11.2. Pessary extractor PRM2 (if required),
 - 11.3. Insert instruction.
12. Ring w / Knob pessary with support consisting of:
 - 12.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7,
 - 12.2. Pessary extractor PRM2 (if required),
 - 12.3. Insert instruction.
13. Hodge's pessary without support, consisting of:

- 13.1. Pessary, size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,
- 13.2. Pessary extractor PRM2 (if required),
- 13.3. Insert instruction.
14. Hodge's pessary with support from:
 - 14.1. Pessary, size 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,
 - 14.2. Pessary extractor PRM2 (if required),
 - 14.3. Insert instruction.
15. Pessary Dish without support consisting of:
 - 15.1. Pessary, size 0, 1, 2, 3, 4, 5, 6, 7,
 - 15.2. Pessary extractor PRM2 (if required),
 - 15.3. Insert instruction.
16. Pessary Dish with support from:
 - 16.1. Pessary size 0, 1, 2, 3, 4, 5, 6, 7,
 - 16.2. Pessary extractor PRM2 (if required),
 - 16.3. Insert instruction.
17. Pessary Cube (Cubes) consisting of:
 - 17.1. Pessary, size 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10,
 - 17.2. Pessary extractor PRM2 (if required),
 - 17.3. Insert instruction.
18. Pessary Cube (Cubes) with drainage consisting of:
 - 18.1. Pessary, size 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10,
 - 18.2. Pessary extractor PRM2 (if required),
 - 18.3. Insert instruction.
19. Pessary Donut (Donuts) consisting of:
 - 19.1. Pessary, size 0, 1, 2, 3, 4, 5, 6, 7,
 - 19.2. Pessary extractor PRM2 (if required),
 - 19.3. Insert instruction.
20. Pessary Donut (Donuts) inflatable consisting of:
 - 20.1. Pessary, size: 2, 3, 4, 5,
 - 20.2. Pessary pear,
 - 20.3. Insert instruction.
21. Pessary Oval (Ovals) without support consisting of:
 - 21.1. Pessary size 1, 2, 3, 4, 5, 6, 7, 8, 9,
 - 21.2. Pessary extractor PRM2 (if required),
 - 21.3. Insert instruction.
22. Ovals pessary with support from:
 - 22.1. Pessary size 1, 2, 3, 4, 5, 6, 7, 8, 9,
 - 22.2. Pessary extractor PRM2 (if required),
 - 22.3. Insert instruction.
23. Pessary Flexi shelf consisting of:
 - 23.1. Pessary, size FS06, FS08, FS10, FS12,
 - 23.2. Pessary extractor PRM2 (if required),
 - 23.3. Insert instruction.
24. Pessary fitting Ring with handle (Ring w / Knob) with FS 8000 support, consisting of:
 - 24.1. Pessary - 6 pcs / unit,
 - 24.2. Pessary extractor PRM2 (if required),
 - 24.3. Insert instruction.
25. Pessary fitting room Maryland (Marland) with FS 7000 support consisting of:
 - 25.1. Pessary - 6 pcs / pack;
 - 25.2. Pessary extractor PRM2 (if required),
 - 25.3. Insert instruction.
26. Pessary fitting Ring (Ring) with FS 6000 support, consisting of:
 - 25.1. Pessary - 6 pcs / unit,

- 25.2. Pessary extractor PRM2 (if required),
- 25.3. Insert instruction.
- 27. Dressing room pessary Gellhorn FS 5000 consisting of:
 - 27.1. Pessary - 6 pcs / unit,
 - 27.2. Pessary extractor PRM2 (if required),
 - 27.3. Insert instruction.
- 28. Fitting pessary Ring FS 1000 consisting of:
 - 28.1. Pessary - 6 pcs / unit,
 - 28.2. Pessary Extractor PRM2 (if required),
 - 28.3. Insert instruction.
- 29. Fitting room pessary Dish with FS 3000 support, consisting of:
 - 29.1. Pessary - 6 pcs / unit,
 - 29.2. Pessary extractor PRM2 (if required),
 - 29.3. Insert instruction.
- 30. Pessary Fitting Cube (Cubes) FS 2000 consisting of:
 - 30.1. Pessary - 6 pcs / unit,
 - 30.2. Pessary extractor PRM2 (if required),
 - 30.3. Insert instruction.
- 31. Fitting pessary Oval with support (Ovals) FS 4000 consisting of:
 - 31.1. Pessary - 6 pcs / unit,
 - 31.2. Pessary extractor PRM2 (if required),
 - 31.3. Insert instruction.

2.1. Manufacturer:

Panpac Medical Corporation, Taiwan, 6F-2, No. 202, Sec.3, Ta-Tong Road, Shi-Chih Dist., New Taipei City, Taiwan, 22103.

Tel.: +886-2-2716-0135, +886-2-2716-0531

e-mail: sales@panpac.com.tw

2.2. Authorized Representative in the Russian Federation

TEKNO-MEDICAL Limited Liability Company, TEKNO-MEDICAL LLC, 119415, Moscow, Vernadsky Prospekt, 39, et. 5 room 1 room 46,

Tel: + 7 (495) 540 5106

e-mail: rostov@tekno-medical.ru


3. Appointment

The urogynecological pessary medical device is intended to restore and maintain the tone of the muscles of the pelvic floor with the aim of treating the pathology of the pelvic organs.

Potential consumer: Professional medical personnel (urologist, gynecologist), as well as patients who have undergone detailed training by the attending physician. The fitting pessaries and the PRM2 pessary extractor are intended for use by a physician only in healthcare facilities.

3.1. Scope

The scope is urology, gynecology.

3.2. Contraindications

- Abnormal uterine bleeding
- Acute genitourinary infections
- Acute pelvic infections
- Patients not complying with the requirements for use.

3.3. Possible side effects

- Bacterial vaginosis,
- Irritation, ulcers in the vagina, possible bleeding,
- Urinary incontinence during exertion (during coughing, sneezing or physical exertion).
- Urinary tract infections
- Difficulties with sexual intercourse (the possibility of injury to a partner).

3.4. Precautionary measures

- Before use, read the instructions for use.
- The doctor must have the skills to install and select the product.
- Each patient should be instructed in detail and trained by a medical specialist in self-handling the product at home (extraction, cleaning, disinfection, administration).
- The product has successfully passed the test for biocompatibility and allergenicity. However, patients with an allergic reaction to silicone should not use this product.
- The doctor must make sure that the patient is not allergic to the material of the product within the first 24-48 hours after installing the pessary. Examine the patient's vagina and find out if the patient had any malaise, chest pain, tenderness, or unusual vaginal discharge. Also determine if there has been an improvement in her personal symptoms.
- The doctor should plan follow-up visits according to the needs of the patient.
- The use of pessaries (with the exception of fitting pessaries) by different patients is not allowed. Pessaries are intended for only one patient (with the exception of the fitting pessaries, the responsibility of the clinic for cleaning, disinfection and sterilization).
- The product is supplied non-sterile. Before use, it must undergo cleaning, disinfection and sterilization in a steam autoclave or chemical sterilization, with the exception of pessaries such as Donuts. (see the section "Cleaning, disinfection, sterilization)."
- Pessaries such as Donuts and Donuts inflatable are not subject to sterilization. Before the first use and each subsequent use, it is subjected to cleaning and disinfection. This type of pessary should be inspected each time by the introduction of cracks.
- Fitting pessaries are designed for different patients. Before the first use and after each fitting, the fitting pessaries should undergo a cleaning, disinfection and sterilization procedure.
- Follow the cleaning, disinfection and sterilization procedures validated by the manufacturer and specified in the instructions for use.
- To facilitate the introduction of the product, the manufacturer recommends treating it with a lubricant-lubricant. (For example, "Gel-lubricant O'KEY® PLUS, tubes of 15 g, 50 g according to TU 9398-001-25074123-13" produced by Biorhythm Laboratory LLC, Russia, Registration certificate No. RZN 2016/4131 dated 23.05.2016).).

4. Performance characteristics, therapeutic indications

Urology-gynecological pessaries are temporary vaginal means of fixation and support of the pelvic organs. They are a medical device of a conservative, non-surgical, treatment method that is inserted into the vagina.

Pessaries are soft, flexible silicone products made from biologically inert silicone of various shapes and sizes. Some types of pessaries to hold the shape are reinforced from the inside with metal inserts. Pessaries without support (without a supporting membrane) restore the correct position of organs by fixing them, but have a lack of probability of squeezing the uterus into the hole of the pessary. Support pessaries have an additional membrane inside the loop. The support membrane not only fixes, but also supports internal organs, preventing the uterus from being squeezed out.

Perforated pessaries (with drainage) provide the best outflow of fluid with increased vaginal secretion.

The product is intended for only one patient (with the exception of the fitting pessaries and the pessary extractor).

The PRM2 pessary extractor is a sleek hook-shaped tool that allows the doctor to easily grab and remove the pessary from the patient's vagina.

To select the appropriate size and shape of the pessary for prolonged use, the doctor should use the fitting pessaries. Before trying on a pessary, the patient empties the bladder.

The fitting pessaries are blue, the pessaries for continuous use are white.

The time of use of the pessary is very individual, it is determined by the attending physician with regular monitoring of the patient. The time of use of the product depends on the speed of restoration of the tone of the muscles of the pelvic floor.

Depending on the shape and size, the product can be shown (see Table 1) with:

- uterine prolapse of 1, 2 and 3 degrees,
- uterine prolapse in combination with stress urinary incontinence,
- cystocele (prolapse of the bladder),
- rectocele (protrusion of the rectum towards the vagina),


Table 1. Therapeutic Indications

Type (option) / with support without support	Prolapse 1 and 2 degrees	Prolapse 2 and 3 degrees	Cystocele	Retrocele	Uterine prolapse combined with stress urinary incontinence
Ring	X	-	X	-	-
Ring w/Knob	X	X	X		X
Ovals	X	-	X	-	-
Shaatz	-	X	X	-	-
Dish	X	X	X	-	
Marland	-	X	X	-	X
Gellhorn with a short, long leg.	-	X	X	X	-
Donuts, Donuts inflatable	-	X	X	X	-
Cubes With drainage / without drainage	-	X	X	X	-
Gehrung	-	X	X	X	-
Hodge	X	-	X	-	X
Cup	X	-	X	-	X
Flexi shelf		X	X	X	-

Table 2.
Coition with the device inside

Possible	Not possible/ before sexual intercourse, you must remove the product
- Ring - Ring w/Knob - Ovals - Shaatz - Dish - Marland - Gehrung - Hodge - Cup	- Gellhorn with a short and long leg, - Donuts - Donuts inflatable, - Cubes with drainage / without drainage, - Flexi shelf

5. Application risks

Product risks are determined in accordance with the requirements of EN ISO 14971: 2012.

As part of the assessment of the risk management analysis, all hazards in relation to risk are acceptable.

In the process of risk analysis, it was concluded that this product has a low level of risk, and any existing risks were eliminated or reduced through safety tests, the right choice of materials and complete instructions for use. It follows that a medical device provides more benefits than harm.

Residual risks are indicated in the Precautions section.


6. Specifications

6.1. Variant of execution of Marland Pessary without support;

size: 2, 3, 4, 5, 6, 7, 8.

The pessary is ring-shaped with a raised platform. The raised platform during installation puts pressure on the neck of the urethra, preventing involuntary urination when coughing and sneezing. It has seven sizes. The pessary is internally reinforced with an elastic element.

Name	Visual appearance	Weights and dimensions					
Marland Pessary without support		Size	Outside diameter ± 1 mm	Inside diameter ± 1 mm	H1 (Height) H2 (Thickness)	Mass ± 2 g	Gross weight ± 3 g
		2	59 mm	35 mm	12 mm 34 mm	32 g	54 g
		3	64 mm	40 mm	12 mm 36 mm	35 g	55 g
		4	70 mm	47 mm	12 mm 38 mm	40 g	60 g
		5	76 mm	52 mm	12 mm 40 mm	45 g	65 g
		6	82 mm	58 mm	12 mm 42 mm	50 g	70 g
		7	88 mm	64 mm	12 mm 44 mm	60 g	80 g
		8	94 mm	70 mm	12 mm 46 mm	68 g	88 g

6.2. Variant of execution Marland Pessary with support;

size: 2, 3, 4, 5, 6, 7, 8.

The pessary is ring-shaped with a raised platform. It has a supporting membrane with drainage holes. It has seven sizes. The pessary is reinforced from the inside by an elastic element.

Name	Visual appearance	Weights and dimensions					
Marland Pessary without support	 The diameter of the four drainage holes 4 ± 0.5 mm	Size	Outside diameter ± 1 mm	Inside diameter ± 1 mm	H1 (Height) H2 (Thickness)	Mass ± 2 g	Gross weight ± 3 g
		2	59 mm	35 mm	12 mm 34 mm	34 g	54 g
		3	64 mm	40 mm	12 mm 36 mm	35 g	55 g
		4	70 mm	47 mm	12 mm 38 mm	40 g	60 g
		5	76 mm	52 mm	12 mm 40 mm	45 g	65 g
		6	82 mm	59 mm	12 mm 42 mm	52 g	72 g
		7	88 mm	65 mm	12 mm 44 mm	60 g	80 g
		8	94 mm	70 mm	12 mm 46 mm	68 g	88 g

6.3. Variant of execution of Shaatz Pessary with drainage;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8.

The pessary has a biconvex disc with drainage holes. It does not have an elastic reinforcing element. It has nine sizes.

Name	Visual appearance	Weights and dimensions					
Shaatz pessary with drainage		Size	Outside diameter ± 1 mm	Inside diameter ± 1 mm	Thickness ± 1 mm	Mass ± 2 g	Gross weight ± 3 g
		0	38 mm	16 mm	10 mm	15 g	35 g
		1	45 mm	16 mm	10 mm	17 g	37 g
		2	51 mm	16 mm	10 mm	20 g	40 g
		3	57 mm	16 mm	10 mm	35 g	45 g
		4	62 mm	16 mm	10 mm	30 g	50 g
		5	70 mm	16 mm	10 mm	40 g	60 g
		6	76 mm	16 mm	10 mm	48 g	69 g
		7	82 mm	16 mm	10 mm	55 g	75 g
		8	88 mm	16 mm	10 mm	65 g	85 g
The diameter of six small holes is 4 ± 0.5 mm (for sizes 0 ~ 5); 5 ± 0.5 mm (for sizes 6 ~ 8)							

6.4. Cup Pessary without support;

size: 0, 1, 2, 3, 4, 5, 6, 7.

Cup Pessary without support has the shape of a bowl, with a large central hole and perforated by six small drainage holes located on the periphery. It does not have a reinforcing element. Presented in eight sizes.


Set Name	Visual appearance	Weights and dimensions					
Cup Pessary without support		Size	Outside diameter ± 1 mm	Inside diameter ± 1 mm	Thickness ± 1 mm	Mass ± 2 g	Gross weight ± 3 g
		0	50 mm	18 mm	18 mm	14 g	34 g
		1	55 mm	23 mm	18 mm	16 g	36 g
		2	57 mm	28 mm	18 mm	18 g	38 g
		3	65 mm	33 mm	18 mm	21 g	41 g
		4	70 mm	38 mm	18 mm	24 g	44 g
		5	75 mm	43 mm	18 mm	25 g	45 g
		6	78 mm	51 mm	18 mm	26 g	46 g
		7	85 mm	52 mm	18 mm	27 g	47 g
		Diameter of six small holes 4 ± 0.5 mm					

6.5. Cup Pessary with support;

size: 0, 1, 2, 3, 4, 5, 6, 7.

The Cup Pessary with support is shaped like a bowl, perforated with six small drainage holes located on the periphery and two small holes in the middle. It does not have a reinforcing element. Presented in eight sizes.

Name	Visual appearance	Weights and dimensions			
Cup Pessary with support	 Diameter of eight small holes 4 ± 0.5 mm (for sizes 0-2, 6-7) 5 ± 0.5 mm (for sizes 3-5)	Size	Outside diameter ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	50 mm	14 g	34 g
		1	55 mm	16 g	36 g
		2	57 mm	18 g	38 g
		3	65 mm	21 g	41 g
		4	68 mm	24 g	44 g
		5	75 mm	25 g	45 g
		6	78 mm	30 g	50 g
		7	85 mm	32 g	52 g

6.6. Gellhorn pessary with a short leg and drainage;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8.

The pessary has the shape of a bowl on a leg. On the base of the leg a thickening is performed. The foot of the pessary is intended for its introduction and retention as deep as possible in the vagina. This pessary model is intended for short-term use - no more than 6-12 hours, administration for a day and removal for a night. Pessary has 9 sizes.


Name	Visual appearance	Weights and dimensions				
Pessary Gellhorn with a short leg and drainage	 Diameter of six small holes 4 ± 0.5 mm Diameter of central hole 6.0 ± 0.5 mm Leg diameter 16 ± 0.5 mm Leg diameter in the wide part 20 ± 1 mm	S i z e	Outside diameter ± 1 mm	H (Height) ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	38 mm	44 mm	20 g	40 g
		1	44 mm	44 mm	24 g	44 g
		2	51 mm	44 mm	27 g	47 g
		3	57 mm	44 mm	35 g	55 g
		4	63 mm	44 mm	40 g	60 g
		5	70 mm	44 mm	50 g	50 g
		6	76 mm	44 mm	55 g	75 g
		7	83 mm	44 mm	60 g	80 g
		8	89 mm	44 mm	65 g	85 g

6.7. Gellhorn pessary with a long leg and drainage;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8.


The pessary has the shape of a bowl on a long leg. On the base of the leg a thickening is performed. The foot of the pessary is intended for its introduction and retention as deep as possible in the vagina. This pessary model is intended for short-term use - no more than 6-12 hours, administration for a day and removal for a night. Pessary has nine sizes.

Name	Visual appearance	Weights and dimensions				
Pessary Gellhorn with long leg and drainage	 Diameter of six small holes 4 ± 0.5 mm (for sizes 0 - 5); 5 ± 0.5 mm (for sizes 6-8) Diameter of center hole 6 ± 0.5 mm Leg diameter 16 ± 0.5 mm Leg diameter in the wide part 20 ± 1 mm	Size	Outside diameter ± 1 mm	Height ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	38 mm	52 mm	22 g	42 g
		1	45 mm	52 mm	26 g	46 g
		2	51 mm	52 mm	28 g	48 g
		3	57 mm	52 mm	36 g	56 g
		4	63 mm	52 mm	40 g	60 g
		5	70 mm	52 mm	52 g	72 g
		6	76 mm	52 mm	56 g	76 g
		7	82 mm	52 mm	62 g	82 g
		8	88 mm	52 mm	67 g	87 g

6.8. Gehrung pessary;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8.

The pessary is made of thin silicone, but reinforced from the inside along the perimeter with an elastic element. Has nine sizes.

Name	Visual appearance	Weights and dimensions				
Gehrung pessary	 Has no holes	Size	Width ± 1 mm	Length ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	45 mm	50 mm	23 g	43 g
		1	50 mm	55 mm	26 g	46 g
		2	50 mm	60 mm	28 g	48 g
		3	50 mm	65 mm	31 g	51 g
		4	55 mm	70 mm	34 g	54 g
		5	55 mm	75 mm	38 g	58 g
		6	60 mm	80 mm	42 g	62 g
		7	60 mm	85 mm	47 g	67 g
		8	65 mm	90 mm	52 g	72 g

**6.9. Ring pessary without support,
size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9.**


The pessary has a ring shape, reinforced with an elastic element inside. The pessary has ten sizes.

Name	Visual appearance	Weights and dimensions					
Ring pessary without support		Size	Outside D ± 1 mm	Inside D ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	44 mm	22 mm	11 mm	14 g	34 g
		1	50 mm	26 mm	12 mm	16 g	36 g
		2	56 mm	32 mm	12 mm	18 g	38 g
		3	63 mm	39 mm	12 mm	20 g	40 g
		4	69 mm	45 mm	12 mm	22 g	42 g
		5	76 mm	52 mm	12 mm	24 g	44 g
		6	83 mm	59 mm	12 mm	28 g	48 g
		7	88 mm	64 mm	12 mm	30 g	50 g
		8	95 mm	65 mm	15 mm	32 g	52 g
		9	100 mm	70 mm	15 mm	34 g	54 g

**6.10. Ring pessary with support,
size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9.**

The pessary has a ring shape with a supporting membrane perforated with six holes. The pessary has ten sizes.


Name	Visual appearance	Weights and dimensions				
Ring pessary with support	 The diameter of the four small drainage holes 4 ± 0.5 mm, Diameter of two large drainage holes 11 ± 0.5 mm Size "0" has only two large holes with a diameter of 12 ± 0.5 mm 	Size	Outside D ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	44 mm	11 mm	16 g	36 g
		1	50 mm	12 mm	17 g	37 g
		2	55 mm	12 mm	18 g	38 g
		3	63 mm	12 mm	22 g	42 g
		4	69 mm	12 mm	24 g	44 g
		5	76 mm	12 mm	27 g	57 g
		6	83 mm	12 mm	30 g	50 g
		7	88 mm	12 mm	35 g	55 g
		8	95 mm	15 mm	37 g	57 g
		9	100 mm	15 mm	42 g	62 g

6.11. Ring pessary with knob (Ring w / Knob) without support

size: : 0, 1, 2, 3, 4, 5, 6, 7,

The pessary is ring-shaped with an oval-shaped thickening and a large central opening. The thickening is intended for additional fixation of the urethra. The pessary is internally reinforced with an elastic element. When installing the pessary, the ring is moved up through the vaginal vault so that the thickening (knob) lifts the jumper between the bladder and the urethra. Has eight sizes


Name	Visual appearance	Weights and dimensions				
Ring pessary with Knob (Ring w / Knob) without support	 H1 22 mm$\pm$ 1 mm H2 29 mm$\pm$ 1 mm	Size	Outside D $\pm$ 1 mm	Thickness $\pm$ 1 mm	Weight $\pm$ 2 g	Gross weight $\pm$ 3 g
		0	44 mm	12 mm	20 g	40 g
		1	50 mm	12 mm	23 g	43 g
		2	58 mm	12 mm	26 g	46 g
		3	63 mm	12 mm	28 g	28 g
		4	69 mm	12 mm	30 g	50 g
		5	76 mm	12 mm	34 g	54 g
		6	83 mm	12 mm	36 g	56 g
		7	88 mm	12 mm	39 g	59 g

6.12. Ring pessary with knob (Ring w / Knob) with support

size: 0, 1, 2, 3, 4, 5, 6, 7,

The pessary has a ring shape with an oval-shaped thickening and a supporting membrane perforated with six holes. The pessary has eight sizes.


Name	Visual appearance	Weights and dimensions				
Ring pessary with Knob (Ring w / Knob) with support	 H1 22 mm $\pm$ 1 mm H2 29 mm $\pm$ 1 mm Diameter of four small drainage holes 4 $\pm$ 0.5 mm Diameter of two large drainage holes 11 $\pm$ 0.5mm Size "0" has only two large drainage holes with a diameter of 12 $\pm$ 0.5mm 	Size	Outside D $\pm$ 1 mm	Thickness $\pm$ 1 mm	Weight $\pm$ 2 g	Gross weight $\pm$ 3 g
		0	44 mm	12 mm	22 g	42 g
		1	50 mm	12 mm	25 g	45 g
		2	57 mm	12 mm	28 g	48 g
		3	63 mm	12 mm	30 g	50 g
		4	69 mm	12 mm	34 g	54 g
		5	76 mm	12 mm	36 g	56 g
		6	83 mm	12 mm	38 g	58 g
		7	88 mm	12 mm	42 g	62 g

6.13. Hodge Pessary without support

size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,


The design of the pessary allows you to give the product any shape. This makes it possible to use the Hodge pessary in situations where the patient's anatomical features do not allow the use of traditional pessaries. Reinforced along the perimeter with an elastic frame. Has 10 sizes

Name	Visual appearance	Weights and dimensions					
Hodge Pessary without support		Size	H1 ± 1 mm	H2 ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	65 mm	38 mm	08 mm	20 g	40 g
		1	70 mm	40 mm	08 mm	22 g	42 g
		2	77 mm	40 mm	08 mm	24 g	44 g
		3	80 mm	45 mm	08 mm	25 g	45 g
		4	85 mm	48 mm	08 mm	26 g	46 g
		5	90 mm	51 mm	08 mm	28 g	48 g
		6	97 mm	53 mm	08 mm	30 g	50 g
		7	100 mm	57 mm	08 mm	35 g	55 g
		8	105 mm	60 mm	08 mm	37 g	57 g
		9	110 mm	63 mm	08 mm	40 g	60 g

6.14. Hodge Pessary with support

Size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9,

The design of the pessary allows you to give the product any shape. This makes it possible to use the Hodge pessary in situations where the patient's anatomical features do not allow the use of traditional pessaries. It has a supporting membrane. Reinforced along the perimeter with an elastic frame. Has 10 sizes.


Name	Visual appearance	Weights and dimensions						
Hodge pessary with support		Size	H1 ± 1 mm	H2 ± 1 mm	H3 ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	65 mm	38 mm	28 mm	08 mm	22 g	42 g
		1	70 mm	40 mm	30 mm	08 mm	24 g	44 g
		2	77 mm	40 mm	32 mm	08 mm	26 g	46 g
		3	80 mm	45 mm	35 mm	08 mm	27 g	47 g
		4	85 mm	48 mm	37 mm	08 mm	28 g	48 g
		5	90 mm	51 mm	38 mm	08 mm	30 g	50 g
		6	97 mm	53 mm	38 mm	08 mm	32 g	52 g
		7	100 mm	57 mm	43 mm	08 mm	37 g	57 g
		8	105 mm	60 mm	45 mm	08 mm	39 g	59 g
		9	110 mm	63 mm	47 mm	08 mm	41 g	61 g

6.15. Dish pessary without support;

size: 0, 1, 2, 3, 4, 5, 6, 7.


The pessary has the shape of a bowl with an oval-shaped thickening and a large central opening. The thickening is intended for additional fixation of the urethra. It has 8 sizes.


Name	Visual appearance	Weights and dimensions						
Dish pessary without support	 Diameter of six small drainage holes $4 \pm 0.5\text{mm}$	Size	Outside D $\pm 1\text{ mm}$	Inside D $\pm 1\text{ mm}$	Thickness $\pm 1\text{ mm}$	H1 $\pm 1\text{ mm}$	Weight $\pm 2\text{ g}$	Gross weight $\pm 3\text{ g}$
		0	50 mm	18 mm	17 mm	10 mm	17 g	27 g
		1	55 mm	23 mm	17 mm	13 mm	19 g	29 g
		2	60 mm	28 mm	17 mm	15 mm	21 g	41 g
		3	65 mm	33 mm	17 mm	17 mm	24 g	44 g
		4	72 mm	38 mm	17 mm	27 mm	25 g	45 g
		5	75 mm	43 mm	17 mm	27 mm	28 g	48 g
		6	80 mm	48 mm	17 mm	27 mm	31 g	51 g
		7	85 mm	52 mm	17 mm	28 mm	32 g	52 g

6.16. Dish pessary with support;

size: 0, 1, 2, 3, 4, 5, 6, 7.

The pessary is in the shape of a bowl with an oval-shaped thickening and a supporting membrane perforated with eight drainage holes. The thickening is intended for additional fixation of the urethra. It has 8 sizes.

Name	Visual appearance	Weights and dimensions					
Dish pessary With support	 Diameter of eight drainage holes $4 \pm 0.5\text{mm}$	Thickness	Outside D $\pm 1\text{ mm}$	Thickness $\pm 1\text{ mm}$	H1 $\pm 1\text{ mm}$	Weight $\pm 2\text{ g}$	Gross weight $\pm 3\text{ g}$
		0	50 mm	17 mm	10 mm	17 g	27 g
		1	55 mm	17 mm	13 mm	19 g	29 g
		2	60 mm	17 mm	15 mm	21 g	41 g
		3	65 mm	17 mm	17 mm	24 g	44 g
		4	72 mm	17 mm	27 mm	25 g	45 g
		5	75 mm	17 mm	27 mm	28 g	48 g
		6	80 mm	17 mm	27 mm	31 g	51 g
		7	85 mm	17 mm	28 mm	32 g	52 g


6.17. Cubes Pessary;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10.

The pessary has a cube shape with concave sides. For ease of removal, the cubic pessaries are equipped with a loop for easier extraction.

This pessary model is intended for short-term use (no more than 6-12 hours, administration for a day and removal for a night). It has eleven sizes.


Name	Visual appearance	Weights and dimensions			
Cubes Pessary	 With loop for easy removal	Size	H1 = H2 ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	25 mm	12 g	32 g
		1	29 mm	14 g	34 g
		2	33 mm	16 g	36 g
		3	37 mm	22 g	42 g
		4	41 mm	28 g	48 g
		5	45 mm	36 g	56 g
		6	50 mm	48 g	72 g
		7	56 mm	60 g	80 g
		8	63 mm	72 g	92 g
		9	70 mm	84 g	104 g
		10	75 mm	96 g	116 g

6.18. Cubes Pessary with drainage;

size: 0, 1, 2, 3, 4, 5, 6, 7, 8, 9, 10.

The pessary has a cube shape with concave sides. For ease of removal, the cubic pessaries are equipped with a loop for easier extraction. The vacuum effect allows you to use this type of product, even if the pelvic floor is not strong enough.

This pessary model is intended for short-term use (no more than 6-12 hours, administration for a day and removal for a night). In the pessary model, silicone cubic perforated perforations are made on the sides for the outflow of vaginal discharge. It has eleven sizes.


Name	Visual appearance	Weights and dimensions			
Cubes Pessary with drainage	 With loop for easy removal Diameter of drainage holes 4 ± 0.5 mm	Size	H1 = H2 ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	25 mm	12 g	32 g
		1	29 mm	14 g	34 g
		2	33 mm	16 g	36 g
		3	37 mm	22 g	42 g
		4	41 mm	28 g	48 g
		5	45 mm	36 g	56 g
		6	50 mm	48 g	72 g
		7	56 mm	60 g	80 g
		8	63 mm	72 g	92 g
		9	70 mm	84 g	104 g
		10	75 mm	96 g	116 g

6.19. Donuts Pessary;

size: 0, 1, 2, 3, 4, 5, 6, 7.


Pessary in the form of a donut filled with air. Not subject to autoclaving. It has eight sizes.


Name	Visual appearance	Weights and dimensions					
Donuts Pessary		Size	Outside D ± 1 mm	Inside D ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		0	50 mm	15 mm	17.5 mm	18 g	38 g
		1	59 mm	16 mm	20.5 mm	22 g	42 g
		2	65 mm	16 mm	23.5 mm	28 g	48 g
		3	69 mm	16 mm	25.5 mm	35 g	55 g
		4	75 mm	16 mm	28 mm	46 g	66 g
		5	83 mm	16 mm	28 mm	56 g	76 g
		6	88 mm	16 mm	35 mm	66 g	86 g
		7	95 mm	18 mm	38.5 mm	75 g	95 g

6.20. Inflatable Donuts pessary;

size 2, 3, 4, 5;

A pessary in the form of a donut is equipped with a pear for filling (pumping) air. Not subject to autoclaving. It has four sizes.


Name	Visual appearance	Weights and dimensions				
Inflatable Donuts pessary	 Squeeze pressure pump length 135 ± 2 mm	Size	Outside D ± 1 mm	Inside D ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		2	51 mm	10 mm	28 g without squeeze pressure pump, 66 g with squeeze pressure pump	82 g
		3	57 mm	11 mm	30 g without squeeze pressure pump, 68 g with squeeze pressure pump	84 g
		4	64 mm	11 mm	32 g without squeeze pressure pump, 70 g with squeeze pressure pump	86 g
		5	70 mm	11 mm	34 g without squeeze pressure pump, 72 g with squeeze pressure pump	88 g

6.21 Ovals Pessary without support;

sizes: 1, 2, 3, 4, 5, 6, 7, 8, 9.


The pessary has an oval shape, reinforced with an elastic element. It has nine sizes.

Name	Visual appearance	Weights and dimensions				
Ovals Pessary without support		Size	H1 ± 1 mm	Thickness s± 1 mm	Weight t ± 2 g	Gross weight ± 3 g
		1	50.5 mm	12 mm	16 g	36 g
		2	57 mm	12 mm	18 g	38 g
		3	63.5 mm	12 mm	20 g	40 g
		4	70.5 mm	12 mm	22 g	42 g
		5	76 mm	12 mm	24 g	44 g
		6	83 mm	12 mm	26 g	46 g
		7	89 mm	12 mm	28 g	48 g
		8	96 mm	14 mm	30 g	50 g
		9	101 mm	14 mm	32 g	52 g

6.22 Ovals Pessary with support;


sizes: 1, 2, 3, 4, 5, 6, 7, 8, 9.

The pessary has an oval shape with a supporting membrane perforated by six drainage holes. Reinforced along the contour with an elastic element. Has nine sizes.


Name	Visual appearance	Weights and dimensions				
Ovals Pessary without support	 The diameter of the four drainage holes is 4 ± 0.5 mm. Two larger holes have a length of 12 ± 1 mm (oval)	Size	H1 ± 1 mm	Thickness s ± 1 mm	Weight ± 2 g	Gross weight ± 3 g
		1	50.5 mm	12 mm	16 g	36 g
		2	57 mm	12 mm	18 g	38 g
		3	63.5 mm	12 mm	20 g	40 g
		4	70.5 mm	12 mm	22 g	42 g
		5	76 mm	12 mm	26 g	46 g
		6	83 mm	12 mm	30 g	50 g
		7	89 mm	12 mm	34 g	54 g
		8	96 mm	14 mm	36 g	56 g
		9	101 mm	14 mm	40 g	60 g


6.23 Flexi shelf pessary with package inserts

The pessary has the shape of a curved bed with a curved foot and two drainage holes.


Name	Visual appearance	Weights and dimensions			
Flexi shelf pessary		Size	H1 ±0.5 mm	Weight ± 2 g	Gross weight ± 3 g
		FS06	57 mm	18 g	36 g
		FS08	64 mm	20 g	38 g
		FS10	70 mm	22 g	40 g
		FS12	76 mm	24 g	42 g
		Diameter of two drainage holes 5 ± 0.5 mm			

6.24. Ring Fitting Pessary with Knob (Ring w / Knob) FS 8000 with support



Name	Visual appearance	Weights and dimensions				
Ring Fitting Pessary with Knob FS 8000 with support	 H1: 23.5 mm$\pm$ 1 mm H2: 28.5 mm$\pm$ 1 mm Diameter of two small drainage holes 4.5 $\pm$ 0.5mm Diameter of two large drainage holes 10 $\pm$ 0.5mm 	Size e	Outside D $\pm$ 1 mm	Thickness s $\pm$ 1 mm	Weight $\pm$ 2 g	Gross weight $\pm$ 10 g
2	58 mm	12 mm	28 g	776 g		
3	63 mm	12 mm	30 g			
4	69 mm	12 mm	34 g			
5	76 mm	12 mm	36 g			
6	83mm	12 mm	38 g			
7	88 mm	12 mm	42 g			

6.25. Marland Fitting Pessary FS 7000 with support

Name	Visual appearance	Weights and dimensions					
Marland Fitting Pessary FS 7000 with support	  Diameter of four holes 4 ± 0.5 mm	Size	Outside D ± 1 mm	Inside D ± 1 mm	Height H1 Thickness H2 ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	58 mm	34 mm	12 mm 34 mm	34 g	852 g
		3	64 mm	40 mm	12 mm 36 mm	35 g	
		4	70 mm	46 mm	12 mm 38 mm	40 g	
		5	76 mm	52 mm	12 mm 40 mm	45 g	
		6	82 mm	58 mm	12 mm 42 mm	52 g	
		7	88 mm	64 mm	12 mm 44 mm	60 g	

6.26. Ring Fitting Pessary FS 6000 with support

Name	Visual appearance	Weights and dimensions			
Ring Fitting Pessary FS 6000 with support	  The diameter of the four small drainage holes is 4.5 ± 0.5 mm, Diameter of two large drainage holes 10 ± 0.5 mm	Size	Outside D ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	58 mm	18 g	714 g
		3	63 mm	22 g	
		4	69 mm	24 g	
		5	76 mm	27 g	
		6	83 mm	30 g	
		7	88 mm	35 g	


6.27. Gellhorn fitting pessary FS 5000


Name	Visual appearance	Weights and dimensions				
Gellhorn fitting pessary FS 5000	 Diameter of six small holes 4 ± 0.5 mm Diameter of the central hole 6 ± 0.5 mm Leg diameter 16 ± 0.5 mm The diameter of the legs in the wide part 20 ± 1 mm 	Size	D ± 1 mm	Height ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	51 mm	52 mm	28 g	1195 G
		3	57 mm	52 mm	36 g	
		4	63 mm	52 mm	40 g	
		5	70 mm	52 mm	50 g	
		6	76 mm	52 mm	56 g	
		7	82 mm	52 mm	68 g	

6.28. Ring Fitting Pessary FS 1000


Name	Visual appearance	Weights and dimensions					
Ring Fitting Pessary FS 1000		Size	Outside D ± 1 mm	Inside D ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	58 mm	32 mm	12 mm	18 g	700 g
		3	63 mm	38 mm	12 mm	20 g	
		4	69 mm	44 mm	12 mm	22 g	
		5	76 mm	50 mm	12 mm	24 g	
		6	83 mm	58 mm	12 mm	28 g	
		7	88 mm	64 mm	12 mm	30 g	



6.29. Dish Fitting Pessary with support for FS 3000

Name	Visual appearance	Weights and dimensions					
Dish Fitting Pessary with support for FS 3000		Size	Outside D ± 1 mm	Thickness ± 1 mm	H1 ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	60 mm	17 mm	15 mm	21 g	778.7 g
		3	65 mm	17 mm	17 mm	24 g	
		4	70 mm	17 mm	27 mm	25 g	
		5	75 mm	17 mm	22 mm	28 g	
		6	80 mm	17 mm	27 mm	31 g	
		7	85 mm	17 mm	28 mm	32 g	
							
Diameter of eight drainage holes 4.5 ± 0.5 mm							

6.30. Cubes fitting Pessary FS 2000

Name	Visual appearance	Weights and dimensions			
Cubes fitting Pessary FS 2000	 	Size	Side length H1 = H2 ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	33 mm	16 g	750 g
		3	37 mm	22 g	
		4	41 mm	28 g	
		5	45 mm	36 g	
		6	50 mm	48 g	
		7	56 mm	60 g	

6.31. Ovals Fitting Pessary with support for FS 4000

Name	Visual appearance	Weights and dimensions				
Ovals Fitting Pessary with support for FS 4000	 The diameter of the four drainage holes is 4 ± 0.5 mm. Two larger holes have a length of 5 ± 1 mm (oval) 	Size	H1 ± 1 mm	Thickness ± 1 mm	Weight ± 2 g	Gross weight ± 10 g
		2	57 mm	12 mm	18 g	720 g
		3	63.5 mm	12 mm	20 g	
		4	70.5 mm	12 mm	22 g	
		5	76 mm	12 mm	26 g	
		6	83 mm	12 mm	30 g	
		7	89 mm	12 mm	34 g	

6.32. Pessary extractor PRM2;

Set Name	Visual appearance	Weights and dimensions		
Pessary extractor PRM2		Size ± 2 mm	Weight ± 2 g	Gross weight ± 3 g
		180 mm	30 g	38 g

6.1. Silicone hardness

Pessaries are made of silicone of varying degrees of hardness.

Type :	Shor hardness ($\pm$ 5A before and after sterilization)
Marland Pessary without support, with support	Frame: 70A Membrane: 40 A,
Shaatz Pessary	40 A
Cup Pessary without support, with support	Frame: 70A
Gellhorn Pessary, short, long-legged	40A
Gehrung Pessary	40A
Ring Pessary without support, with support, Ring Pessary (Ring w / Knob) with Knob	Frame, Knob: 70A Membrane: 40 A,
Hodge pessary	40A
Dish Pessary, without support, with support	Frame, thickening knob: 70A
Cubes pessary	30A
Donuts Pessary, Inflatable Donuts pessary	50A
Oval pessary with support, without support	Frame 70A Membrane 40 A,
Flexi shelf Pessary	70A

7. Accessories Description

No accessories.

8. Drug content

The product does not contain any medicinal products, or materials of animal or human origin.


9. Mode of application.

9.1. Commissioning

- The product is delivered in individual packaging, non-sterile.
- Before the first use, the product is removed from the packaging and subjected to cleaning, disinfection and sterilization by a doctor in a healthcare facility by autoclaving or by chemical methods (with the exception of pessaries such as Donuts and Donuts inflatable), which are not suitable for sterilization and are only cleaned and disinfected)
- Immediately before the introduction of the pessary must be treated with a small amount of lubricant-lubricant (for example, O'KEY® PLUS Lubricant Gel, tubes of 15 g, 50 g according to TU 9398-001-25074123-13 "manufactured by Biorhythm Laboratory, Russia, Registration Certificate No. RZN 2016/4131 dated 23.05.2016).

The sterilized and grease-treated product is ready for use.


Pessaries Donuts and Donuts inflatable are ready for use after disinfection and treatment with a lubricating lubricant.

9.2. Cleaning, disinfection, sterilization before first use

The manufacturer validated two sterilization methods for the product: steam and chemical, confirmed by the document "Evaluation of the effectiveness of sterilization of the product urogynecological pessary by chemical or steam method."

Name	Method of cleaning, disinfection, sterilization, sterilization before first use
Pessaries (not fitting ones, for individual use)	
Marland without support, Marland with support, Shaatz with drainage, Cup without support, Cup with support Gellhorn with a short leg and drainage, Gellhorn with a long leg and drainage, Gehrung, Ring without support Ring with support Ring with a knob (Ring w / Knob) without support, Ring with a knob (Ring w / Knob) with support, Hodge without support, Hodge with support Dish without support Dish with support, Cube Cubes with drainage, Ovals without support Ovals pessary with support, Flexi shelf	<u>In the health care facility by medical personnel:</u> I Method Cleaning: Neutral detergent, soft brush, warm running water. Disinfection: Alcohol-containing disinfectants - 75% ethanol solution (place the product for 20 minutes). Steam Sterilization - Type B Autoclave, Temperature 121°C (do not exceed 124°C) Sterilization exposure for 20 minutes, Overpressure 0.21 ± 0.01 MPa Products withstand 150 cycles without depreciation of quality. II method Cleaning: Neutral detergent, soft brush, warm running water. Disinfection, sterilization: CIDEX OPA disinfectant (SIDEX OPA) manufactured by Advance Sterilization Products (USA) Johnson & Johnson Corporation (USA). (RU.77.99.88.002.E.002947.04.13 04/15/2013) 20-minute exposure time
Fitting pessaries	
Ring with a knob (Ring w / Knob) with FS 8000 support Marland FS 7000 Ring FS 6000 with support Gellhorn FS 5000 Ring FS 1000 Dish FS 3000 Cubes FS 2000 Oval FS 4000	<u>in health care facilities</u> I. method Cleaning: Neutral detergent, soft brush, warm running water. Disinfection: Alcohol-containing disinfectants - 75% ethanol solution (place the product for 20 minutes). Sterilization: 1. Type B steam sterilizer, Temperature 121°C (do not exceed 124°C) Sterilization exposure for 20 minutes, Overpressure 0.21 ± 0.01 MPa Products withstand 150 cycles without depreciation of quality. II method. Cleaning: Neutral detergent, soft brush, warm running water. Disinfection, chemical sterilization: CIDEX OPA disinfectant manufactured by Advanced Sterilization Products (USA) of Johnson & Johnson Corporation (USA).

	(RU.77.99.88.002.E.002947.04.13 04/15/2013) 20-minute exposure time
Pessary type Donuts and Donuts inflatable	
Donuts and Inflatable Donuts pessary	Not to be sterilized, only disinfection Cleaning: Neutral detergent, soft brush, warm running water. Disinfection: 75% ethanol solution (place the product for 20 minutes).
Squeeze pressure pump for Inflatable Donuts pessary	Not to be cleaned, disinfected or sterilized.
Pessary Extractor PRM2	
Pessary Extractor PRM2	<u>in health care facilities</u> Cleaning: Neutral detergent, soft brush, warm running water. Disinfection: Alcohol-containing disinfectants - 75% ethanol solution (place the product for 20 minutes). Sterilization: Type B steam autoclave Temperature 121°C (do not exceed 124°C) Sterilization exposure for 20 minutes, Overpressure 0.21 ± 0.01 MPa Products withstand 150 cycles without depreciation of quality.

9.3. Selection of product size

Attention! The procedure for selecting the type and size of the pessary is carried out only by a professional medical professional. The doctor must have the ability to administer pessaries.

Attention! Before trying on and introducing the pessary the doctor must conduct a study of the pelvic floor of the patient, make sure that there are no contraindications for the use of pessaries. Pelvic examination also helps determine the appropriate type and size of pessary.

Attention! Before trying on the product, the patient must empty the bladder.


To begin the fitting, the doctor must evaluate the width of the vagina and use this information to select a pessary of suitable size. The patient should be recommended the largest pessary, which is suitable for her, but does not cause inconvenience. The patient is examined in a prone position and standing, with a sample of Valsalva and without it.

The examining physician should be able to position the finger on both sides of the pessary. Choosing the size of the pessary, the doctor focuses on the fact that during installation he does not squeeze the bladder and intestines, but does not fall out during movement, coughing and stress.

After several days or weeks, the doctor may adjust the selection of the product. After the correct type and size of the pessary is selected and successfully installed, the patient must independently remove and reinstall the pessary


9.4. Pessary installation

Attention! Before administration of the pessary, the prolapsed organs must be returned to the anatomically correct position.

1. To determine the size of the pessary, the medical specialist uses the fitting pessaries. The pessary should be optimal, not too large and not too small so as not to fall out or squeeze the cervix excessively.

2. Before the introduction of the pessary, the doctor compresses the product (Fig. A) and in a compressed form inserts it into the patient's vagina (Fig. B). The pessary is inserted as deep as possible into the vagina (Fig. B). Thanks to the elastic properties of silicone, the pessary straightens and takes its original shape.


3. To check the correct size and location of the pessary, the doctor must manually examine the distance between the outer edge of the pessary and the vaginal wall. The interval in the thickness of the finger will provide patient comfort and reduce the risk of necrosis from pessary pressure. With the right size of the pessary, a small gap should remain between the product and the vaginal wall (10-15 mm, finger width).

4. If the interval is insufficient, the following smaller size should be used.

5. If the interval between the vaginal wall and the pessary is too large (more than the thickness of the doctor's finger), the pessary may be too small. In this case, it is recommended to use the next larger pessary.

6. The sequence of introduction of the pessary Donuts inflatable is shown in the table below:

Fig.1	Connect the valve of the pessary hose (B) to the connector (C) of the blue bulb (A) and completely drain the air.
Fig. 2	The pessary should take the form shown in Figure 2. Disconnect hose valve (B) from bulb valve (C). An air valve installed in the pessary's hose will not allow the pessary to take its original shape.
Fig 3	4. Carefully insert the squeezed pessary into the vagina. The hose should protrude beyond the vagina.

 Fig.4	5. Connect the pessary hose (B) to the plug (D) of the white bulb. Inflate the pessary by pressing the pear (A), Attention! Do not pump the pessary; it may malfunction and cause injury to the patient.
 Fig.5	6. Disconnect bulb (A) from valve (B). Insert the hose into the vaginal vault. Use your thumb and middle finger to grip the pessary, and then carefully remove the pessary by slowly pulling it out of the vagina. Attention! DO NOT pull on the hose when removing the pessary. This can damage the pessary. For personal hygiene, the pessary should be removed and cleaned daily.


7. After installation, any type of pessary should not move when standing, sitting on a chair, squatting, tilting down. The patient should not feel discomfort with any routine movements. Correct installation should be assessed with tension in the abdominal muscles (e.g., coughing, straining) and with movement.

8. No additional tools are required when installing the pessary. The doctor must manually enter the product using sterile gloves. With proper installation, the patient does not feel discomfort or pain. When extracting the pessary, the doctor can use the pessary extractor. Pessaries Cubes are equipped with hinges for ease of extraction.


9. After removing the fitting pessary, it is necessary to clean, disinfect, sterilize it.

10. To prevent the occurrence of possible side effects after installing the pessary, the patient must undergo a gynecological examination with a cytological examination of the discharge once every three months or more often.


The schematic arrangement of various types of pessaries in the vagina is shown in the figure below:

A) For supporting pessaries such as Ring, Ovals, Marland, Cup, Gehrung, Hodge.

B) For pessaries of the types Cubes, Donuts, Gellhorn, Flexi shelf, Dish.

9.5. Female Patient Instruction

Attention! The patient should be instructed in detail and trained by a medical specialist in self-care with the pessary at home (extraction, administration, daily cleaning, disinfection).


Attention! Patients daily (or 2 times a day) should carry out cleaning and disinfection of the used pessary home. (see section 7.2. "Cleaning, disinfection, sterilization").

Attention! The patient should regularly inspect the Donuts pessary and Donat (Donuts) inflatable for the presence / absence of cracks or scratches on the surface. A damaged pessary is not to be used! If there are cracks or scratches, the pessary must be disposed of in accordance with "Disposal").

For self-installation of the pessary, a woman can put one foot on a chair or the edge of the bed and in this position insert the pessary into the vagina. If it is too difficult, you can enter the pessary in a prone position.

10. Instructions for reprocessing and sterilization

Attention	- Cleaning, disinfection, sterilization of the fitting pessaries and the pessary extractor is carried out by a trained medical specialist. The specialist should be in a work coat, headgear, rubber gloves. - All pessaries selected for use by the same patient, when re-processed, are cleaned and disinfected at home by the patient herself, according to the instructions for use. - The patient should remove and carry out cleaning and disinfection of the pessary daily (or 2 times a day). - For cleaning, disinfection and sterilization, use only the products and methods specified in the instructions for use.
Restrictions on reprocessing	Retreatment has minimal impact on the product. The fitting pessaries and the pessary extractor can maintain 150 sterilization cycles without depreciation of quality.
INSTRUCTIONS	
Place of use	No special requirements
Protection and transportation	No special requirements
Training	No special requirements. The product is not collapsible
Manual cleaning	- Rinse the product thoroughly with a neutral detergent and a brush under warm running water. - Dry the product with a clean cloth.
Disinfection	To disinfect all pessaries and the pessary extractor with an alcohol-containing disinfectant - 75% ethanol solution (place the product for 20 minutes).
Drying	1-3 minutes to complete evaporation of the alcohol-containing agent
Packaging	The product is placed in a suitable size craft bag for autoclaving. For the chemical sterilization method, packaging is not required.
Sterilization	Only fitting pessaries and pessary extractor are re-sterilized. I method Type B steam autoclave Temperature 121°C (do not exceed 124°C) Sterilization exposure for 20 minutes, Overpressure 0.21 ± 0.01 MPa Products withstand 150 cycles without depreciation of quality. II method. Chemical sterilization CIDEX OPA disinfectant (SIDEX OPA) manufactured by Advance Sterilization Products (USA) Johnson & Johnson Corporation (USA). (RU.77.99.88.002.E.002947.04.13 04/15/2013) 20-minute exposure time
Storage	- To determine the sterilization validity after autoclaving in the autoclave, refer to the instructions for use of the manufacturer of the used kraft bag. - After chemical sterilization, it is necessary to store the product in a sterile sterilization box lined with a sterile sheet for a period of not more than 3 days. - After disinfection at home, the product is placed in a clean plastic bag or container until next use.
Additional	Secondary treatment (cleaning, disinfection and sterilization) of the fitting pessaries and the pessary

Information	extractor is carried out by a medical specialist in a hospital. Secondary treatment (cleaning, disinfection) of pessaries selected by a doctor for one patient is carried out by the patient at home.
Manufacturer	Panpac Medical Corporation, Taiwan, 6F-2, No. 202, Sec.3, Ta-Tong Road, Shi-Chih Dist., New Taipei City, Taiwan, 22103.

11. Pessary extraction

- **Attention!** If the procedure for extracting the pessary is carried out independently by the patient at home, you must thoroughly wash your hands before this procedure.
- Cubes pessary for ease of extraction is equipped with a silicone loop, which is available for finger grip.
- Pessaries without a loop should be gripped with your fingers, squeezed slightly, and squeezed out of the vagina. Due to the elastic properties of the silicone, the discomfort when removing the pessary is minimal.
- For self-extraction of the pessary, the patient can put one foot on a chair or the edge of the bed. If it is too difficult, you can enter the pessary in a prone position.
- **Attention!** The PRM2 pessary extractor may only be used by a physician.

12. Duration of pessaries use

The duration of treatment with the help of a pessary is determined by the doctor and depends on a combination of many factors (stability of the vaginal microbiocenosis, degree of estrogen deficiency) and, therefore, is very individual. On average, the patient can use the pessary for several months, subject to constant monitoring by the attending physician.

13. Circumstances requiring medical advice

The patient should immediately contact her doctor in the following cases:

- Occurrence of discomfort / pain when wearing a pessary
- Appearance of discharge of yellow or greenish color, indicating the presence of pathogenic bacteria
- Itching sensation in the vagina
- Difficulties with the introduction / removal of the product.

14. Terms of Use

The medical device is intended for use at temperatures from +10 to + 43 ° C, air humidity from 10 to 100%. It can be used in hospitals and at home.

15. Information contained in the leaflet

- Each product is accompanied by a brief package insert containing the following information:
- Appointment
- Indications
- Contraindications
- Possible side effects
- Instructions for reprocessing and sterilization
- Circumstances requiring medical advice

16. Material Information

Product name or parts and accessories (researched)	Materials used in the manufacture
1. Marland Pessary without support 2. Marland Pessary with support 3. Shaatz pessary with drainage 4. Cup Pessary without support 5. Cup Pessary with support. 6. Gellhorn pessary with a short leg and drainage 7. Gellhorn pessary with a long leg and drainage 8. Gehrung pessary 9. Ring pessary without support. 10. Ring pessary with support 11. Ring pessary with a knob (Ring w / Knob) without support. 12. Ring pessary with knob (Ring w / Knob) with support. 13. Hodge pessary without support. 14. Hodge pessary with support. 15. Dish pessary without support. 16. Dish pessary with support 17. Cube Pessary Cube. 18. Cubes pessary with drainage. 19. Donuts pessary 20. Inflatable Donuts pessary 21. Ovals pessary without support 22. Ovals pessary with support. 23. Flexi shelf Pessary	Silicone - brand XIAMETER, manufacturer Dow Corning Co, Silicone pigment # 1924, white manufacturer Shenzhen Kanglibang Science and Technology Co., Ltd
24. Ring Fitting Pessary with a knob (Ring w / Knob) FS 8000, 25. Marland fitting pessary with FS 7000 support 26. Ring Fitting Pessary FS 6000 with support 27. Gellhorn Fitting Pessary FS 5000 28. Ring fitting Pessary FS 1000 29. Dish Fitting pessary FS 3000 30. Cubes Fitting pessary FS 2000 31. Ovals Fitting pessary with support FS 4000	Silicone brand "XIAMETER" manufactured by Dow Corning Co, Silicone pigment # 1526, Cobalt Blue, manufacturer Shenzhen Kanglibang Science and Technology Co., Ltd,
32. Squeeze pressure pump for Inflatable Donuts pessary	Polyvinyl chloride SG-3, Cromofix black ebony dye
33. Pessary Extractor PRM2	Stainless steel Jis G 4303, Coating: Zibo Nylon Blue, manufactured by Yusheng Plastic Products Co. Ltd.

17. Environmental requirements

The product during transportation, storage and use does not have a negative impact on the environment.

18. Applicable international standards

ISO 13485:2016,

MDD 93/42/EEC as amended by 2007/47/EC,

EN ISO 14971:2012,

ISO 15223-1:2012.

19. Transportation

Products are transported in a cardboard box by all means of transport in covered vehicles with protection from direct sunlight, in accordance with the rules of transport applicable to this type of transport at an air temperature of -10 °C to + 60 °C and air humidity from 10 to 100%.

20. Storage conditions, shelf life, service life

Storage Environment:

air temperature from + 0 °C to + 60 °C,

humidity from 10 to 100%,

Protect from direct sunlight.

Shelf life of a medical device:

Shelf life is 5 years from the date of packaging and is indicated on the packaging for the following options:

- Extractor pessaries PRM2,
- Inflatable Donuts pessary,
- Flexi shelf Pessary,
- all fitting pessaries.

Shelf life is 10 years from the date of packaging and is indicated on the packaging for the following options:

Marland without support, Marland with support, Shaatz with drainage, Cup without support, Cup with support, Gellhorn with short leg and drainage, Gellhorn with long leg and drainage, Gehrung, Ring without support, Ring with support, Ring with knob (Ring w / Knob) without support, Ring with knob (Ring w / Knob) with support, Hodge without support, Hodge with support, Dish without support, Dish with support, Cube, Cube with drainage, Donuts, Ovals without support, Ovals pessary with support.

The life of a medical device is determined by the number of sterilization cycles. The product can withstand without compromising its properties more than 150 sterilization cycles (autoclaving or chemical sterilization).

21. Packaging and labeling of a medical device

21.1. Individual packing

- Pessaries for one patient and are packaged in individual packaging made of cardboard and shrink PVC film, in which a brief package insert is enclosed. The size of an individual package of pessaries for one patient: Width 13.5 ± 0.3 cm, Length 12.5 ± 0.3 cm. Height depends on the size of the product, but does not exceed 8.0 cm.
- Pessary Donuts inflatable is packed in a plastic bag with a ZIP-LOCK clasp measuring 32 ± 0.5 cm x 14 ± 0.5 cm.
- The fitting pessaries are packed in transparent plastic containers, separated by partitions for easy transportation.
 - The size of the plastic container for the Gellhorn fitting pessary FS 5000 is $(22.0 \times 19.0 \times 7.0) \pm 0.3$ cm,
 - The size of the plastic container for all other types of fitting pessaries is: $(15.0 \times 10.5 \times 10.5) \pm 0.3$ cm.
- Plastic containers are packed in boxes of three-layer corrugated cardboard.
 - Size $(25.5 \times 21.5 \times 8.5) \pm 0.5$ cm for the Gellhorn fitting pessary FS 5000,
 - Size: $(16.5 \times 12.0 \times 12.5) \pm 0.5$ cm for all other types of fitting pessaries.

- The pessary extractor is packed in a transparent bag made of polypropylene with an adhesive valve, the size of the closed bag $(20.0 \times 10.0) \pm 0.5$ cm.

Each individual pessary package contains a package insert. The weight of the instruction insert is 6 ± 1 g.

Individual packaging of a pessary is a consumer container. The presence or absence of a pessary extractor (supplied if necessary) is marked on the pessary packaging with a "V" (presence) or "-" (absence) sign.

Each individual package contains an insert leaflet. The mass of the package insert is 6 ± 1 g.

21.2. Individual packaging labeling

This marking applies to all versions of a medical device. It contains:

- Name of product and version
- Product size (if available)
- Serial number
- packing date
- Date "Use until"
- Name and address of official distributor
- Name and address of manufacturer
- Name of authorized representative in the EU
- Product Conformity Mark to EU Directive
- Symbol "Refer to the instructions for use"
- Symbol "Not sterile"

21.3. Individual packaging labeling project for the Russian Federation

Name of MD:

RU № XXXXXXXX from XX.XX.XXXX

Design variant: XXXXXXXX

Size (if any)

Pessary extractor PRM2

X.XXX.X

00.00.0000

00/00/0000

Name, address

2460

Symbols on the label

	Batch code
	Date of manufacture
	Use up
	The product is not sterile
	Caution! Refer to the instructions for use.
	Refer to the instructions for use.
	Manufacturer (with address).
	CE mark
	Keep out of sunlight
	Transportation temperature limit
	Mark on the presence / absence of the PRM2 Pessary Extractor in the product

21.4. Labeling layout for the Russian Federation


Photo labels for RF

22. Disposal

Used products are disposed of as class B medical waste in accordance with the rules of SanPiN 2.1.3684-21 (subject to mandatory decontamination (neutralization), disinfection).

Unused products after the expiration of the storage period are disposed of:
as class A medical waste.

23. Maintenance and Repair Requirements

Products are not subject to maintenance and repair.

24. Warranty

The warranty shelf life is 5 years from the date of packaging and is indicated on the packaging for the following options:

- Extractor pessaries PRM2,
- Inflatable Donuts pessary,
- Flexi shelf pessary
- all fitting pessaries.

The warranty shelf life is 10 years from the date of packaging and is indicated on the packaging for the following options:

Marland without support, Marland with support, Shaatz with drainage, Cup without support, Cup with support, Gellhorn with short leg and drainage, Gellhorn with long leg and drainage, Gehrung, Ring without support, Ring with support, Ring with Knob (Ring w / Knob) without support, Ring with knob (Ring w / Knob) with support, Hodge (Hodge) without support, Hodge with support,

Dish without support, Dish with support, Cube, Cube with drainage, Donuts, Ovals without support, Ovals pessary with support.

The warranty period of storage begins from the moment of delivery of the goods is the residual period from the date of packaging.

Goods with broken packaging and goods with identified manufacturing defects should be presented for replacement within 30 days from the delivery of the selling organization

The buyer has the right to make complaints regarding the quality of the goods. Customer satisfaction extends to medical devices that are used only for their intended purpose and in accordance with the instructions for use. Improper use of the medical device is not provided. The manufacturer does not bear any responsibility for any possible damage caused to the consumer if the latter does not comply with the requirements of the manufacturer specified in the instructions for use.

25. Complaints

Claims on the quality of goods within the territory of the Russian Federation can be sent by the buyer to the authorized representative in the Russian Federation LLC TEKNO-MEDICAL, 119415, Moscow, Vernadsky Prospekt, 39, 5th floor, prem. 5, off. 1,

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