
Medtronic, Inc.
710 Medtronic Parkway N.E.,
Minneapolis MN 55432,
USA.
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879


Declaration Letter


Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, Technical Manual for SelectSecure™ MRI SureScan™ model 3830 is true and unaltered copies of the original documents.

Name: Melanie Vincent
Title: Senior Regulatory Affairs Director
NOTARY STAMP

28 Oct 2021
Date

28 October 2021

Medtronic

SelectSecure™ MRI SureScan™ 3830

Steroid eluting, bipolar, implantable, nonretractable screw-in, atrial/ventricular, catheter delivered, transvenous lead

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The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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Medtronic, SelectSecure, SureScan

Contents

1	Device description	3
2	Indications for use	3
3	Contraindications	4
4	Warnings and precautions	4
5	Potential complications	5
6	Directions for use	6
7	Detailed device description	10
8	CE mark of conformity	11
9	Explanation of symbols on package labeling	11
10	Medtronic disclaimer of warranty	12
11	Service	12

1 Device description

The Medtronic SelectSecure MRI SureScan Model 3830 steroid eluting, bipolar, implantable, nonretractable screw-in, atrial/ventricular, catheter delivered, transvenous lead is designed for pacing and sensing in the atrium or ventricle.

The lead has been tested for use in the Magnetic Resonance Imaging (MRI) environment. The following lead lengths are MR conditional: 59 cm, 69 cm, and 74 cm. Other lead lengths are not MR conditional.

The lead has a nonretractable helical electrode made of titanium nitride coated platinum alloy for active fixation in the endocardium by rotating the lead body in a clockwise direction. Active fixation leads are particularly beneficial for patients who have smooth or hypertrophic hearts where lead dislodgement may be a potential problem.

The lead also has a second, larger electrode made of titanium nitride coated platinum alloy proximal to the tip electrode and an IS-1¹ Bipolar (BI) connector. The lead features MP35N nickel alloy conductors, silicone inner insulation, and polyurethane outer insulation.

The distal tip contains a maximum of 100.0 µg of beclomethasone dipropionate. Upon exposure to body fluids, the steroid elutes from the lead tip. Steroid is known to suppress the inflammatory response that is believed to cause threshold rises typically associated with implanted pacing electrodes.

Note: To implant the Model 3830, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and an introducer valve which allows passage through or removal from an IS-1 connector. Contact your local Medtronic representative for further information regarding compatible delivery systems.

¹ IS-1 BI refers to an International Connector Standard (ISO 5841-3) whereby pulse generators and leads so designated are assured of a basic mechanical fit.

1.1 Medtronic SureScan system

A complete SureScan system is required for use in the MR environment. A complete SureScan system includes a Medtronic SureScan device with the appropriate number of Medtronic SureScan leads. Any other combination may result in a hazard to the patient during an MRI scan.

The Model 3830 lead is part of the Medtronic SureScan system. The SureScan system includes a Medtronic SureScan device connected to Medtronic SureScan leads. Labeling for SureScan system components displays the SureScan logo and the MR Conditional symbol.

SureScan logo

MR Conditional symbol. The Medtronic SureScan system is MR Conditional and is designed to allow implanted patients to undergo an MRI scan under the specified MRI conditions for use.

The MRI SureScan feature permits a mode of operation that allows a patient with a SureScan device to be safely scanned by an MRI machine while the device continues to provide appropriate pacing. When programmed to On, MRI SureScan operation disables arrhythmia detection, magnet mode, and all user-defined diagnostics. **Before performing an MRI scan, refer to the SureScan system MRI technical manual for important information about procedures and MRI-specific warnings and precautions.**

1.2 Package contents

The lead and accessories are supplied sterile. Each package contains the following items:

- 1 lead with anchoring sleeve
- 1 vein lifter
- Product documentation

1.3 Accessory descriptions

Dispose of all single-use accessories according to local environmental requirements.

Anchoring sleeve – An anchoring sleeve secures the lead from moving and protects the lead insulation and conductors from damage caused by tight sutures.

Vein lifter – A vein lifter facilitates catheter or introducer insertion into a vessel.

2 Indications for use

The Model 3830 lead has application where implantable atrial or ventricular, single-chamber or dual-chamber pacing systems are indicated. The Model 3830 lead is intended for pacing and sensing in the atrium or ventricle.

English 3

3 Contraindications

The following are contraindications for use of Medtronic implantable, screw-in, transvenous leads.

- Use of ventricular transvenous leads is contraindicated in patients with tricuspid valvular disease.
- Use of ventricular transvenous leads is contraindicated in patients with mechanical tricuspid heart valves.
- Use of steroid eluting transvenous leads is contraindicated in patients for whom a single dose of 100.0 µg beclomethasone dipropionate may be contraindicated.

4 Warnings and precautions

A complete SureScan system is required for use in the MRI environment. Before performing an MRI scan, refer to the SureScan technical manual for MRI-specific warnings and precautions.

Note: Medical procedure warnings and precautions that pertain to the Medtronic implanted system are provided in the manual that is packaged with the device or on the Medtronic Manual Library website (www.medtronic.com/manuals).

Inspecting the sterile package – Inspect the sterile package with care before opening.

- Contact your local Medtronic representative if the seal or package is damaged.
- Do not store this product above 40°C (104°F).
- Do not use the product after its expiration date.

Sterilization – Medtronic has sterilized the package contents with ethylene oxide before shipment. This lead is for single use only and is not intended to be resterilized.

Single use – The lead and accessories are for single use only.

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during acute lead system testing, the implant procedure, or whenever arrhythmias are possible or intentionally induced during post-implant testing.

Line-powered and battery-powered equipment – An implanted lead forms a direct current path to the myocardium. During lead implant and testing, use only battery-powered equipment or line-powered equipment specifically designed for this purpose to protect against fibrillation that may be caused by alternating currents. Line-powered equipment used in the vicinity of the patient must be properly grounded. Lead connector pins must be insulated from any leakage currents that may arise from line-powered equipment.

Concurrent devices – Output pulses, especially from unipolar devices, may adversely affect device sensing capabilities. If a patient requires a separate stimulation device, either permanent or temporary, allow enough space between the leads of the separate systems to avoid interference in the sensing capabilities of the devices. Previously implanted pulse generators and implantable cardioverter defibrillators should generally be explanted.

Steroid use – It has not been determined whether the warnings, precautions, or complications usually associated with injectable beclomethasone dipropionate apply to the use of this highly localized, controlled-release lead.

Pregnancy – Beclomethasone dipropionate (BDP) has been shown to be teratogenic in many species when given in doses equivalent to the human dose. There are no adequate and well-controlled studies in pregnant women. BDP should be used during pregnancy only if the potential benefit justifies the potential risk to the fetus. Studies in mice, rats, and rabbits have shown that adrenocorticoids increase the incidence of cleft palate, placental insufficiency, and spontaneous abortions, and can decrease the intrauterine growth rate.

Nursing mothers – Systemically administered corticosteroids appear in human milk and could suppress growth, interfere with endogenous corticosteroid production, or cause other untoward effects in nursing infants. Because of the potential for serious adverse reactions in nursing infants from corticosteroids, a decision should be made whether to discontinue nursing or to use a non-steroidal lead, taking into account the importance of the lead and the drug to the mother.

Handling the steroid tip – Avoid reducing the amount of steroid available prior to lead implant. Reducing the available amount of steroid may adversely affect low-threshold performance.

- Do not allow the electrode surface to come in contact with surface contaminants.
- Do not wipe or immerse the electrode in fluid, except blood, at the time of implant.

Handling the lead – Handle the lead with care at all times.

- Do not implant the lead if it is damaged. Return the lead to your Medtronic representative.
- Do not attempt to straighten or realign the helix if the helix is deformed. Return the lead to your Medtronic representative.
- Protect the lead from materials that shed small particles such as lint and dust. Lead insulators attract these particles.
- Handle the lead with sterile surgical gloves that have been rinsed in sterile water or a comparable substance.
- Do not severely bend, kink, or stretch the lead.
- Do not apply pressure to the helix.
- Do not immerse the lead in mineral oil, silicone oil, or any other liquid, except blood, at the time of implant.
- Do not use surgical instruments to grasp the lead.
- Do not force the lead if resistance is encountered during lead passage.
- Keep the helix within the guide catheter of the delivery system if passing through the tricuspid valve to prevent damage to the helix, valve, and/or endocardial tissue.

Vessel and tissue damage – Use care when positioning the lead.

- Consider another site for lead placement other than the apex of the right ventricle if there is reason to believe the patient has an unusually thin wall at the apex of the right ventricle.
- Avoid known infarcted or thin ventricular wall areas to minimize the occurrence of perforation and dissection.

- Avoid acute trauma to the endocardium, including possible perforation, caused by excessive torque and/or tip pressure. Acute trauma to the endocardium may result in temporarily high impedance or threshold values.

Chronic repositioning or removal – Proceed with extreme caution if a lead must be removed or repositioned. Chronic repositioning or removal of screw-in transvenous leads may not be possible because the helix may become deformed and/or entangled as a result of manipulating the lead. In most clinical situations, it is preferable to abandon unused leads in place. Return all removed or unused leads, or lead sections to Medtronic for analysis.

- Observe the helix via fluoroscopy or x-ray before attempting to reposition to determine if the helix shape is intact. If the helix appears deformed, removal may be difficult and is not recommended.
- Lead removal may result in avulsion of the endocardium, valve, or vein.
- Lead junctions may separate, leaving the lead tip and bare wire in the heart or vein.
- Chronic repositioning of a lead may adversely affect a steroid lead's low-threshold performance.
- Abandoned leads should be capped to avoid transmitting electrical signals.
- Severed leads should have the remaining lead end sealed and the lead body sutured to adjacent tissue.

Chronic lead removal and the SureScan pacing system – When implanting a SureScan pacing system, consider the risks associated with removing previously implanted leads before doing so. Abandoned leads or previously implanted non-SureScan labeled leads compromise the ability to safely scan the SureScan pacing system during MRI scans.

Magnetic resonance imaging (MRI) – An MRI is a type of medical imaging that uses magnetic fields to create an internal view of the body. If certain criteria are met and the warnings and precautions provided by Medtronic are followed, patients with an MR Conditional device and lead system are able to undergo an MRI scan; for details, refer to the SureScan MRI technical manual that Medtronic provides for an MR Conditional device.

Diathermy treatment (including therapeutic ultrasound) – Diathermy is a treatment that involves the therapeutic heating of body tissues. Diathermy treatments include high frequency, short wave, microwave, and therapeutic ultrasound. Except for therapeutic ultrasound, do not use diathermy treatments on cardiac device patients. Diathermy treatments may result in serious injury or damage to an implanted device and lead system. Therapeutic ultrasound (including physiotherapy, high intensity therapeutic ultrasound, and high intensity focused ultrasound), is the use of ultrasound at higher energies than diagnostic ultrasound to bring heat or agitation into the body. Therapeutic ultrasound is acceptable if treatment is performed with a minimum separation distance of 15 cm between the applicator and the

implanted device and lead system, as long as the ultrasonic beam is pointing away from the device and lead system.

5 Potential complications

The potential complications related to the use of transvenous leads include, but are not limited to, the following patient-related conditions that can occur when the lead is being inserted and/or repositioned.

- Cardiac perforation
- Cardiac tamponade
- Fibrillation and other arrhythmias
- Heart wall or vein wall rupture
- Infection
- Muscle or nerve stimulation
- Pericardial rub
- Pneumothorax
- Thrombotic and air embolism
- Thrombosis
- Valve damage (particularly in fragile hearts)

Other potential complications related to the lead and the programmed parameters include, but are not limited to, the complications listed in the following table. Symptoms of the following complications include loss of capture or intermittent or continuous loss of capture or sensing²:

Complication	Corrective action to be considered
Lead dislodgement	Reposition the lead
Lead conductor or helix electrode fracture or insulation failure	Replace the lead. In some cases with a bipolar lead, the implantable device may be programmed to a unipolar configuration or the lead may be unipolarized.
Threshold elevation or exit block	Adjust the implantable device output. Replace or reposition the lead

Potential acute/chronic complications associated with lead placement that may require lead replacement to correct include, but are not limited to, the following:

Implant technique	Potential complication
Forcing the lead through the guide catheter	Helix electrode damage and/or insulation damage
Use of too medial an approach with the guide catheter resulting in clavicle & first rib binding	Conductor coil fracture, insulation damage
Puncturing the periosteum and/or tendon when using subclavian guide catheter approach	Conductor coil fracture, insulation damage

² Transient loss of capture or sensing may occur for a short time following surgery until lead stabilization takes place. If stabilization does not occur, lead dislodgement may be suspected.

6 Directions for use

Warning: Before implanting a SureScan system, consider the risks associated with removing previously implanted leads. Abandoned leads or previously implanted non-SureScan labeled leads compromise the ability to safely scan the SureScan system during MRI scans.

Note: To implant the Model 3830, a compatible delivery system is required, such as a Medtronic delivery system. A compatible delivery system includes a guide catheter and an introducer valve which allows passage through or removal from an IS-1 connector. Contact your local Medtronic representative for further information regarding compatible delivery systems.

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. Some implant techniques vary according to physician preference and the patient's anatomy or physical condition.

6.1 Preparing the delivery system

Prepare the delivery system for lead implant according to the instructions in the product documentation packaged with the delivery system.

6.2 Selecting an insertion site


Caution: When using a subclavian approach for insertion, use a more lateral approach to minimize the risk of first rib clavicular crush. First rib clavicular crush may subsequently fracture the lead body.

Caution: Certain anatomical abnormalities, such as thoracic outlet syndrome, may pinch and subsequently fracture the lead body.

The guide catheter assembly may be inserted through several different venous routes, including the right or left cephalic vein, other subclavian branches, or the external or internal jugular vein.

Select an insertion site. See (Figure 1) for suggested insertion site.

Figure 1.

1 Suggested entry site

6.3 Inserting the guide catheter assembly

The guide catheter assembly may be inserted using either the vein lifter in the lead package or the method described in the delivery system product documentation.

The guide catheter assembly may be inserted using venotomy through several venous routes, including the right or left cephalic vein, other subclavian branches, or the external or internal jugular vein.

Insert the guide catheter assembly using the vein lifter:

1. Insert the tapered end of the vein lifter into the incised vein (Figure 2).

Figure 2.

2. Gently push the tip of the guide catheter assembly underneath the vein lifter and into the vein.

See the delivery system product documentation for the recommended method of inserting the guide catheter assembly.

6.4 Positioning the guide catheter

See the delivery system product documentation for details about positioning the guide catheter in the right atrium or right ventricle.

6.5 Inserting the lead into the guide catheter

Warning: For leads that will be placed in the right ventricle, keep the helix within the guide catheter when passing through the tricuspid valve to prevent damage to the helix, valve, and/or endocardial tissue.

Caution: If wiping the lead is necessary prior to insertion, avoid snagging the helix in gauze and ensure that the anchoring sleeve remains in position.

Insert the lead into the guide catheter. Pass the lead through the introducer valve to minimize the backflow of blood.

Note: Keep the helix within the catheter's distal tip.

6.6 Positioning the lead in the ventricle

Warning: To minimize the occurrence of perforation and dissection, avoid known infarcted or thin ventricular wall areas.

Warning: If there is reason to believe the patient has an unusually thin wall at the apex of the right ventricle, the implanter may wish to consider another site for placement of the lead.

Warning: Excessive torque and/or tip pressure may cause acute trauma to the endocardium, including possible perforation. The acute trauma may result in temporary high impedance or threshold values.

Position the lead in the ventricle:

1. Position the tip of the guide catheter in the ventricle. See the delivery system product documentation for details about positioning the guide catheter in the ventricle.
2. When the distal tip of the guide catheter is near the desired location for lead placement, gently advance the lead while retracting the guide catheter until the anode ring is extended approximately 2-3 cm beyond the distal opening of the guide catheter.
3. Use fluoroscopy to facilitate accurate lead placement. Use fluoroscopy (lateral position) to ensure that the tip is not in a retrograde position or is not lodged in the coronary sinus.
4. Rotate the lead body clockwise to affix the helix in the endocardium. It is recommended that the implanter turn the lead body approximately 3 to 4 complete (360°) rotations so that the helix is fully imbedded in the endocardium.
Note: Use either the lead serial number label or the anchoring sleeve to visually count the number of turns while using fluoroscopy.

Accurate positioning of the electrode is essential for stable pacing and sensing. A satisfactory position is usually achieved when the above procedures are followed.

6.7 Positioning the lead in the atrium

Warning: Excessive torque and/or tip pressure may cause acute trauma to the endocardium, including possible perforation. The acute trauma may result in temporarily high impedance or threshold values.

Position the lead in the atrium:

1. Position the tip of the guide catheter in the atrium. See the delivery system product documentation for details about positioning the guide catheter in the atrium.
2. When the distal tip of the guide catheter is near the desired location for lead placement, gently push the lead through the guide catheter until the fixation mechanism is outside the distal opening of the guide catheter.
3. Use fluoroscopy to facilitate accurate lead placement. Generally, a satisfactory lead position has the lead tip situated against the atrial endocardium in or near the apex of the appendage. As viewed on fluoroscopy (A-P view), the lead tip points medially and forward toward the left atrium. A successful position is usually achieved with an anterior, medial, or lateral tip location.

4. Rotate the lead body clockwise to affix the helix in the endocardium. It is recommended that the implanter turn the lead body approximately 4 to 5 complete (360°) rotations so that the helix is fully imbedded in the endocardium.

Note: Use either the lead serial number label or the anchoring sleeve to visually count the number of turns while using fluoroscopy.

If properly positioned, the lead tip will sway from side to side with each atrial contraction (viewed under A-P fluoroscopy). In the absence of spontaneous atrial activity, movement can be produced by pacing the atrium through the lead.

6.8 Verifying helix electrode fixation

Verify helix electrode fixation:

1. Verify fixation using the proper technique depending on whether the lead was fixated in the ventricle or atrium:
 - a. **For a ventricular lead:** Gently pull back on the lead and check for resistance to verify fixation. A properly fixated helix will remain in position. If the helix is not properly fixated, the lead tip will move into the right atrium or may become loose.
 - b. **For an atrial lead:** Gently pull back on the lead and check for resistance to verify fixation. Use frontal fluoroscopy to check for lateral "to-and-fro" movement of the atrial tip, which reflects atrial and ventricular contractions. Poor fixation is suspected when lead tip movement seems random.
2. After confirmation of helix fixation, lead slack is built up in the atrium or ventricle to prevent tip dislodgement. Enough slack is assumed present, if, under fluoroscopy, an atrial lead assumes an "L" shape and a ventricular lead assumes an "S" shape during deep inhalation. Avoid excessive slack buildup that may cause the loop of the atrial lead to drop near the tricuspid valve.
3. Obtain electrical measurements to verify satisfactory placement and electrode fixation. Refer to Section 6.9.
4. If a lead must be repositioned or removed, proceed with caution. Refer to Section 6.10.

6.9 Taking electrical measurements

Take electrical measurements:

1. Pull the guide catheter out of the chamber where the lead is positioned so the guide catheter does not interfere with electrical measurements.
2. Attach the clips of a surgical cable to the lead connector pin and connector ring.
3. Use an implant support instrument to obtain electrical measurements. Medtronic recommends using a pacing system analyzer. For information on the use of the implant support instrument, see the product documentation for that device.

Satisfactory lead placement is indicated by low stimulation thresholds and adequate sensing of intracardiac signal amplitudes.

- A low stimulation threshold provides for a desirable safety margin, allowing for a possible rise in thresholds that may occur within 2 months following implant.
- Adequate sensing amplitudes ensure that the lead is properly sensing intrinsic cardiac signals. Minimum signal requirements depend on the implantable device's sensitivity capabilities. Acceptable acute signal amplitudes for the lead must be greater than the minimum device sensing capabilities, including an adequate safety margin to account for lead maturity.

Table 1. Recommended measurements at implant

Measurement required	Ventricle	Atrium
Maximum acute stimulation thresholds ^a	1.0 V 3.0 mA	1.5 V 4.5 mA
Minimum acute sensing amplitudes	5.0 mV	2.0 mV

^a At a pulse duration setting of 0.5 ms.

4. If electrical measurements do not stabilize to acceptable levels, repositioning and repeating the testing procedure may be necessary. Refer to "Repositioning or removing the lead."

Note: Initial electrical measurements may deviate from the recommendations because of acute cellular trauma. If this occurs, wait 5 to 15 minutes and repeat the testing procedure. Values may vary depending upon lead type, device settings, cardiac tissue condition, and drug interactions.

6.9.1 Checking diaphragmatic stimulation

Diaphragmatic stimulation should also be checked by pacing at 10 V and observing on fluoroscopy whether the diaphragm contracts with each paced stimulus. If diaphragmatic pacing occurs, reduce the voltage until a diaphragmatic pacing threshold is determined. If the diaphragmatic threshold is less than the required programmed pacing output, the lead should be repositioned. Refer to "Repositioning or removing the lead."

6.9.2 Taking pacing impedance (resistance) measurements

Pacing impedance (or resistance) is used to assess device function and lead integrity during routine device patient follow-up and to assist in troubleshooting suspected lead failures. Additional troubleshooting procedures include ECG analysis, visual inspection, measurement of thresholds, and electrogram characteristics.

Pacing impedance values are affected by many factors including lead position, electrode size, conductor design and integrity, insulation integrity, and the patient's electrolyte balance. Apparent pacing impedance is also significantly affected by the measurement technique. Comparison of pacing impedance should be done using consistent methods of measurements and equipment.

An impedance higher or lower than the typical values is not necessarily a conclusive indication of a lead failure. Other causes must be considered as well. Before reaching a conclusive diagnosis, the full clinical picture must be considered. The full clinical picture includes pacing artifact size and morphology changes in 12-lead analog ECGs, muscle stimulation with bipolar leads, sensing and/or capture problems, patient symptoms, and device characteristics.

Recommendations for clinically monitoring and evaluating leads in terms of impedance characteristics are listed below.

For implantable devices with telemetry readout of impedance:

- Routinely monitor and record impedance values at implant and follow-ups using consistent output settings.
Note: Impedance values may be different at different programmable output settings (e.g., pulse width or pulse amplitude) of the implantable device or pacing system analyzer.
- Establish a baseline chronic impedance value once the impedance has stabilized, generally within 6 to 12 months after implant.
- Monitor for significant impedance changes and abnormal values.
- Where impedance abnormalities occur, closely monitor the patient for indications of pacing and sensing problems. The output settings used for measuring impedance should be the same as that used for the original measurements.
- For patients at high risk, such as implantable device-dependent patients, physicians may want to consider further action such as increased frequency of monitoring, provocative maneuvers, and ambulatory ECG monitoring.

For implantable devices without telemetry:

- Record impedance value at implant. Also record the measurement device, its output settings, and the procedure used.
- At the time of device replacement, if pacing system analyzer system-measured impedance is abnormal, carefully evaluate lead integrity (including thresholds and physical appearance) and patient condition before electing to reuse the lead.
- Impedances below 250 Ω may result in excessive battery current drain, which may seriously compromise device longevity, regardless of lead integrity.

For more information on obtaining electrical measurements, consult the product documentation supplied with the testing device.

6.10 Repositioning or removing the lead

Warning: Removal of the lead may result in avulsion of the endocardium, valve, or vein. In addition, the lead junctions may separate leaving the helix or a portion of the lead in the heart or vein. In most clinical situations, it is preferable to abandon unused leads in place. Repositioning the lead after it has been implanted may adversely affect a steroid lead's low-threshold performance.

Caution: If you determine that the lead requires repositioning, consider the possibility that the helix may become deformed and/or entangled as a result of manipulating the lead. Observe the helix via fluoroscopy or x-ray before attempting to reposition to determine if the helix shape is intact. If the helix appears deformed, removal may be difficult and is not recommended.

Reposition or remove the lead:

1. Rotate the lead body counterclockwise to withdraw the helix from the implant site if the helix appears intact and repositioning is required.

Note: The number of counterclockwise rotations needed to withdraw the helix from the implant site before applying traction may be greater than the number of revolutions required for fixation.

Note: If the helix is still imbedded in the endocardium, additional turns on the lead body should be applied rather than applying a retraction force.

2. Counterclockwise rotation should be continued throughout the repositioning and removal process to decrease the possibility of damage to the cardiovascular tissue.
3. Repeat the ventricular or atrial positioning procedure and the verifying helix electrode fixation procedure. Use the guide catheter to reposition the lead.

6.11 Removing the guide catheter from the lead

Once the lead is in the final position, remove the guide catheter from the lead before surgical closure. See the delivery system product documentation for details. Repeat electrical measurements; see Section 6.9.

6.12 Anchoring the lead

Caution: Use care when anchoring the lead.

- Use an anchoring sleeve with all leads.
- Do not use absorbable sutures to anchor the lead.
- Do not secure the sutures so tightly that they damage the vein, lead, or anchoring sleeve.
- Do not use the anchoring sleeve tabs for suturing (Figure 3).
- Do not tie a suture directly to the lead body (Figure 4).
- Do not dislodge the lead tip.
- Do not attempt to remove or cut the anchoring sleeve.
- Do not remove the tabs on anchoring sleeves. Tabs are provided to minimize the possibility of the sleeve entering the vein.
- Do not allow passage of the anchoring sleeve into the guide catheter and/or the venous system.

Figure 3.

1 Tab


Figure 4.

Anchor the lead using all 3 grooves:

1. Position the anchoring sleeve against or near the vein.
2. Secure the anchoring sleeve to the lead body by tying a suture firmly in each of the 3 grooves (Figure 5).

Figure 5.

3. Use at least 1 additional suture in 1 of the grooves to secure the anchoring sleeve and lead body to the fascia.

6.13 Connecting the lead

Connect the lead to the implantable device according to the instructions in the device manual.

Connect the lead to the device.


1. Obtain final electrical measurements.
2. Insert the lead connector into the connector block on the device. For instructions on proper lead connections, see the product documentation packaged with the device.

6.14 Placing the device and lead(s) into the pocket

Caution: Use care when placing the device and lead(s) into the pocket.

- Ensure that the lead(s) does not leave the device at an acute angle.
- Do not grip the lead(s) or device with surgical instruments.
- Do not coil the lead(s) (Figure 6). Coiling the lead(s) can twist the lead body and may result in lead dislodgment.

Figure 6.

Caution: To prevent undesirable twisting of the lead body, wrap the excess lead length loosely under the device and place both into the subcutaneous pocket.

Place the device and lead(s) into the pocket:

1. Rotate the device to loosely wrap the excess lead length under the device (Figure 7).

Figure 7.

2. Insert the device and lead(s) into the pocket.
3. Suture the pocket closed.

6.15 Post-implant evaluation

After implant, monitor the patient's electrocardiogram until the patient is discharged. If a lead dislodges, it typically occurs shortly after implant.

Recommendations for verifying proper lead positioning include x-rays and pacing and sensing threshold measurements.

In the event of a patient death, explant all implanted leads and devices and return them to Medtronic with a completed Product Information Report form. Call the appropriate phone number on the back cover if there are any questions on product handling procedures.

7 Detailed device description


7.1 Specifications (nominal)

Parameter	Model 3830
Type	Bipolar
Chamber	Atrium/Ventricle
Fixation	Nonretractable screw-in
Length	59 cm, 69 cm, 74 cm
Connector	IS-1 BI
Material	Conductors: MP35N
	Connector pin: Stainless steel
	Connector ring: Stainless steel
	Inner insulator: Silicone rubber/ETFE
	Outer insulator: Polyurethane
Electrode material	Helix: Titanium nitride coated platinum alloy
	Ring: Titanium nitride coated platinum alloy
Electrode surface area	Helix: 3.6 mm ²
	Ring: 16.9 mm ²
Tip to ring spacing	9.0 mm

Parameter	Model 3830
Lead body diameter	1.4 mm (4.1 Fr)
Catheter introduction size	3.0 mm (9.0 Fr)
Helix length (exposed)	1.8 mm
Resistance	Unipolar: 29 ± 6 Ω (69 cm)
	Bipolar: 99 ± 22 Ω (69 cm)
Steroid	Beclomethasone dipropionate
Amount of steroid	100.0 µg max

7.2 Specifications drawing (nominal)

Figure 8.

- 1 Helix electrode surface area: 3.6 mm²
- 2 Tip to ring spacing: 9 mm
- 3 Ring electrode surface area: 16.9 mm²
- 4 Anchoring sleeve
- 5 Lead length: 59 cm, 69 cm, 74 cm
- 6 IS-1 BI connector

8 CE mark of conformity

The following CE mark of conformity applies to the Model 3830 lead.

CE0123
2003

9 Explanation of symbols on package labeling

Refer to the package labels to see which symbols apply to this product.

Table 2. Explanation of symbols on package labeling

Symbol	Explanation
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union acts.
	Do not use if package is damaged
	Do not reuse
	Upper limit of temperature
	Sterilized with ethylene oxide
	Consult instructions for use at this website
	Date of manufacture
	Manufacturer
	Use by
	Authorized representative in the European community
	Reorder number
	Serial number
	Lead length

Table 2. Explanation of symbols on package labeling (continued)





Symbol	Explanation
	Open here
	Package contents
	Product documentation
	Accessories
	Inner diameter
	Outer diameter
	Lead
	Transvenous ventricular lead
	Transvenous atrial lead
	Pace
	Sense
	Steroid eluting
	Non-retractable screw-in
	Non-retractable screw-in, Steroid eluting BDP
	Lead introducer with guide wire
	Catheter
	Catheter delivered

Table 2. Explanation of symbols on package labeling (continued)

Symbol	Explanation
	MR Conditional
	SureScan logo

10 Medtronic disclaimer of warranty

For complete disclaimer of warranty information, see the accompanying disclaimer of warranty document.

11 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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
fax 001-763-514-4879

Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, Technical Manual for Deflectable catheter system SelectSite™ C304-S59, C304-L69, C304-XL74 is true and unaltered copies of the original documents.


Name: Melanie Vincent

Title: Senior Regulatory Affairs Director

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English

1 Device description

The device features a percutaneous needle and syringe to access the venous insertion site, a guide wire to access the vein, a valve to reduce blood loss during the implant procedure, a deflectable catheter to introduce a transvenous device, a catheter dilator to facilitate deflectable catheter passage, and a guide catheter splitter to remove the deflectable catheter.

The deflectable catheter features a deflecting distal section, controlled by the deflectable catheter handle, for placement of transvenous devices. The body of the deflectable catheter is radiopaque for visibility on fluoroscopy.

1.1 Contents of package

Each package contains the following items:

- 1 deflectable catheter (either the 6227DEF, C304-S59, C304-L69 or C304-XL74)
- 1 catheter dilator
- 1 valve
- 1 guide catheter splitter
- 1 guide wire with tip straightener
- 1 percutaneous needle
- 1 syringe
- product documentation

2 Indications for use

The device is indicated to provide a pathway through which diagnostic and therapeutic transvenous devices are introduced within the chambers and coronary vasculature of the heart and for introducing balloon catheters into the coronary sinus or leads into vessels of the left heart via the coronary sinus.

3 Contraindications

Use of the device is contraindicated for patients with obstructed or inadequate vasculature and for right ventricular use in patients with tricuspid valvular disease or a mechanical tricuspid heart valve.

4 Warnings and precautions

Transvenous device compatibility – Use the deflectable catheter only with compatible transvenous devices. No test data is available to demonstrate compatibility of the deflectable catheter with any device not manufactured by Medtronic.

- For the Attain 6227DEF Catheter, no test data is available to demonstrate compatibility of transvenous devices with an outer diameter larger than 2.1 mm (8.2 Fr), with the exception of the 2.3 mm (7.0 Fr) catheter dilator included in this package.
- For the C304 Catheters, no test data is available to demonstrate compatibility of the transvenous devices with an outer diameter larger than 1.5mm (4.8Fr), with the exception of the 1.85mm (5.6Fr) dilator included with this package.

Consequences of using the deflectable catheter with incompatible devices may include the inability to deliver the transvenous device or damage to the transvenous device during delivery.

Single use – The device is for single use only and is not intended to be re-sterilized.

Thrombus formation – Thrombogenicity evaluations were conducted using a heparinized model. If your patient cannot be adequately anticoagulated, it is unknown whether thrombus formation may occur with this product.

Inspecting the sterile package – Inspect the sterile package with care before opening it. Contact a Medtronic representative if the seal or package is damaged.

Handling the deflectable catheter – Handle the deflectable catheter with care at all times.

- Do not kink, stretch, or severely bend the deflectable catheter.
- Do not use surgical instruments to grasp the deflectable catheter.
- Avoid contact with liquids other than isopropyl alcohol, blood, saline, or contrast solution.

Using the valve – Use the valve to impede the back flow of venous blood during the implant procedure.

- Ensure that the flush port stopcock is closed before attaching the valve to the deflectable catheter hub.

Necessary hospital equipment – Keep external defibrillation and backup pacing readily available during implant. Use of the device, transvenous devices, or both may cause heart block.

5 Potential complications

Potential complications related to the use of the deflectable catheter may include, but are not limited to, the following events:

- Air embolism
- Allergic reaction to contrast media
- Arteriovenous fistula formation
- Bleeding at the insertion site
- Brachial plexus injury
- Cardiac tamponade
- Dislodgement
- Dissection
- Endocarditis
- Heart block
- Hematoma formation
- Hemothorax
- Infection
- Irregular heart beat
- Mediastinal widening
- Perforation
- Pneumothorax
- Subclavian artery puncture
- Thrombophlebitis
- Thrombosis
- Valve damage
- Vascular occlusion
- Vessel damage

6 Directions for use

Proper surgical procedures and sterile techniques are the responsibility of the medical professional. The following procedures are provided for information only. Each physician must apply the information in these instructions according to professional medical training and experience.

6.1 Preparing the deflectable catheter for insertion

Warning: The deflectable catheter must be thoroughly flushed with saline and free of air prior to use to avoid an air embolism in the patient.

Warning: Do not use a power injection syringe to inject contrast solution.

1. Remove the components of the sterile package within the sterile field.
2. Verify that the deflectable catheter deflects properly by rotating the handle clockwise (Figure 1). Each time the handle is rotated, the distal section is further deflected. After the rotating has stopped, the distal section will hold the amount of deflection for the number of times the handle was rotated.


Figure 1. Deflecting and straightening the deflectable catheter

3. Verify that the deflectable catheter straightens properly by rotating the handle counterclockwise.

Note: After the deflectable catheter has been deflected and straightened, the deflectable catheter may not be completely straight. If necessary, insert a catheter dilator and carefully manipulate the distal section of the deflectable catheter to straighten it as much as possible. The catheter dilator prevents the deflectable catheter shape from becoming deformed.

4. Attach the valve to the hub of the deflectable catheter (Figure 2).


Figure 2. Attaching the valve to the deflectable catheter

5. Flush the deflectable catheter through the valve flush port with saline or heparinized saline prior to use.
6. Flush the catheter dilator.
7. Insert the catheter dilator through the valve and into the deflectable catheter. The distal end of the catheter dilator should extend past the distal end of the deflectable catheter.
8. Set the deflectable catheter assembly aside within the sterile field.

6.2 Inserting the deflectable catheter assembly

Warning: After the deflectable catheter assembly has been advanced into the right atrium, remove or retract the catheter dilator to minimize the risk of perforation.

Caution: To avoid damaging tissue and the tip of the deflectable catheter assembly, keep the catheter dilator fully inserted into the deflectable catheter during passage to the right atrium.

Caution: If the subclavian approach is required for deflectable catheter assembly insertion, avoid techniques that may damage the deflectable catheter assembly.

- The insertion site should be as far lateral as possible to avoid clamping the deflectable catheter assembly between the clavicle and the first rib.
- Do not force the deflectable catheter assembly if significant resistance is encountered during insertion or passage.
- Do not use techniques such as adjusting the patient's posture to facilitate passage. If resistance is encountered, it is recommended to use an alternate venous entry site.

The deflectable catheter assembly may be inserted through several different venous routes, including the right or left cephalic vein or the external or internal jugular vein. Use the cephalic vein whenever possible to avoid damage to the deflectable catheter assembly in the first rib/clavicular (thoracic inlet) space.

Insert the deflectable catheter assembly:

1. Attach the needle to the syringe.
2. Obtain venous access using the needle and syringe.
3. Remove the syringe while holding the needle position fixed when venous blood is aspirated.
4. Straighten the J-shaped tip of the guide wire by advancing the tip straightener over the end of the guide wire.
5. Insert the tapered end of the tip straightener into the needle and feed the guide wire through the needle and into the vein.
6. Remove the tip straightener and advance the guide wire until the distal tip reaches the right atrium. If the guide wire must be withdrawn while the needle is still inserted, both the wire and needle should be withdrawn as a unit.
7. Remove the needle.
8. Place the deflectable catheter assembly over the proximal end of the guide wire.
9. Advance the deflectable catheter assembly through the insertion vein and then to the right atrium.
10. Remove the catheter dilator, leaving the guide wire and deflectable catheter in place.

6.3 Positioning the deflectable catheter in the right atrium or right ventricle

1. If positioning the deflectable catheter in the right ventricle, pass the guide wire through the tricuspid valve. Then use the guide wire to pass the deflectable catheter through the tricuspid valve.
2. Position the deflectable catheter toward the desired location in the right atrium or right ventricle.
3. Use fluoroscopy to verify that the deflectable catheter is in the desired position.

6.4 Positioning the deflectable catheter in the coronary sinus

1. Rotate the deflectable catheter posteriorly and to the patient's left to enter the coronary sinus ostium. Use fluoroscopy and inject small amounts of contrast solution to verify that the deflectable catheter is in the desired position.
2. Advance the deflectable catheter into the coronary sinus.

6.5 Obtaining a venogram

When placing a lead in the coronary sinus, venograms should be obtained to indicate the anatomy of the coronary vasculature.

- When using the Attain 6227DEF catheter, refer to the product documentation packaged with a 2.3 mm (7.0 Fr) or smaller venogram balloon for further instructions on obtaining a venogram.
- When using the SelectSite C304-S59 / C304-L69 / C304-XL74 catheter, refer to the product documentation packaged with a 1.85 mm (5.6 Fr) or smaller venogram balloon for further instructions on obtaining a venogram.

6.6 Inserting the transvenous device into the deflectable catheter and positioning the transvenous device

Insert the transvenous device into the deflectable catheter. After obtaining the desired amount of deflection for placement of the transvenous device, refer to the product documentation packaged with the transvenous device for further positioning instructions. If necessary, use the deflectable catheter to reposition the transvenous device.

6.7 Removing the deflectable catheter from the transvenous device

Warning: Before attempting to remove the deflectable catheter, straighten the distal section as much as possible to avoid vessel damage during removal.

Warning: Keep the cutting surface of the slitter in line with the hub and catheter shaft axis. Do not move the catheter or slitter in a back and forth motion. This ensures that the hub and guide catheter wall are cleanly slit

by the slitter blade, preventing transvenous device damage or dislodgements.

Caution: During and after slitting the deflectable catheter, handle the slit edges of the deflectable catheter carefully to prevent transvenous device damage.

Note: Remove the deflectable catheter with care. Removing the deflectable catheter may dislodge the transvenous device. If the transvenous device dislodges, refer to the product documentation supplied with the transvenous device for repositioning instructions.

When the transvenous device is in final position, the deflectable catheter must be removed from the transvenous device before surgical closure. The deflectable catheter is removed from the transvenous device using a slitter.

Remove the deflectable catheter from the transvenous device:

1. Remove the valve from the deflectable catheter hub (Figure 3).


Figure 3. Removing the valve

2. Follow the appropriate directions as provided in Appendix A for information on slitting.
3. Discard the deflectable catheter and slitter after the deflectable catheter has been removed from the body of the transvenous device.

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§ Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate address or telephone number listed on the back cover.

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



Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, Instructions for use for delivery catheter C315H20, C315H40, C315HIS, C315J, C315S4, C315S5, C315S10 is true and unaltered copies of the original documents.

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USA

Caution: Federal Law (USA) restricts this device to sale by or on the order of a physician.

CE
0050
2010

M715407B001.Jc.sm 7/12/10 4:52 pm
UC2005smrca EN US40US EE MR
4.875 x 7.125 inches (123.825 mm x 180.875 mm)

Medtronic Confidential
Acc_DVD_R04

N715407B001A.book Page 1 Monday, July 12, 2010 4:50 PM

M715407B001 Rev A

Refer to the "Acc_DVD" category in
doc#163256 for Pinning Instructions.

The following are trademarks of Medtronic:

Medtronic

Explanation of symbols on package labeling

Refer to the package label to see which symbols apply to this product

Conformité Européenne (European Conformity). This symbol means that the device fully complies with European Directive AIMD 90/385/EEC.

Date of manufacture

Manufacturer

Use by

Lot number

Sterilized using ethylene oxide

Do not reuse

Maximum storage temperature

Open here/Peel here

Catalogue number

Do not use if package is damaged

For US audiences only

Consult instructions for use

Outer diameter

Manufactured in

EC authorized representative

Inner diameter

Package contents

Guide-catheter dilator

Product documentation

Contents of Package

This package contains product documentation, 1 vessel dilator, and 1 of the following C315 catheters:

C315 H20 & H40

C315 H15

C315 J

C315 S4

C315 S5

C315 S10

M715407B001_sym.dwg 7/12/10 4:53 pm
UC2005xxxxx EN US/OS EE MR
4.875 x 7.125 inches (123.825 mm x 180.875 mm)

Electronic Confidential
Acc_DVD_R04

M715407B001A.book Page 4 Monday, July 12, 2010 4:52 PM

M715407B001 Rev A

Refer to the "Acc_DVD" category in
doc#183256 for Printing Instructions.

1 Device description

The C315 is a single-use percutaneous catheter intended to introduce various types of pacing or defibrillator leads and catheters. The C315 is packaged with a vessel dilator for introduction into the vasculature. Proximally, the C315 is equipped with a hemostatic valve, and the distal tip is radiopaque to facilitate imaging under fluoroscopy. The C315 is designed to be slidable, thereby allowing its removal after device placement. A variety of curves are available to accommodate various anatomies and different lead locations.

1.1 Contents of package

See page 3 for package contents.

2 Indications for use

The C315 is indicated for the introduction of various types of pacing or defibrillator leads and catheters.

3 Contraindications

Use of the C315 percutaneous catheter is contraindicated in patients with obstructed or inadequate vasculature.

4 Warnings and precautions

Necessary hospital equipment – Keep external defibrillation equipment nearby for immediate use during insertion, placement, acute lead system testing, or whenever arrhythmias are possible or intentionally induced. Backup pacing should be readily available during implant. Use of a delivery system and/or leads may cause heart block. For further information about this and other potential adverse events or complications, refer to the technical manual packaged with the appropriate product.

Vessel and tissue damage – Use care when passing a catheter through vessels and tissue.

- Avoid damage, such as perforations and dissections, to vessels and cardiac tissue, during catheter passage and positioning.
- Do not push, pull, or rotate the catheter against resistance. If resistance is met, discontinue movement, determine the reason for resistance, and take appropriate action before continuing. Damage to the catheter may prevent it from performing with accurate torque response and control, and may cause vessel damage.

Device compatibility – Use the catheter only with compatible transvenous devices. No test data is available to demonstrate compatibility of any non-Medtronic device with the catheter. Consequences of using the catheter with incompatible devices may include the inability to deliver the transvenous device or damage to the transvenous device or catheter during delivery.

Slitting evaluations on this product were conducted using Medtronic Universal II or Medtronic Adjustable Slitters. These are the only slitters recommended for use with the C315 catheter.

Inspecting the sterile packages – Inspect the packages prior to opening.

- Contact your local Medtronic representative and provide the lot number from the outer package label if the seal or package is damaged or if there are any concerns regarding product integrity.
- Do not use the product after its expiration date.
- The recommended storage temperature for this product is at or below 40 °C (104 °F).
- The catheter set has been sterilized with ethylene oxide prior to shipment.

For single use only – This device is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize this device for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Handling the catheter – Handle the catheter with care at all times.

- Do not kink, stretch, or severely bend the catheter.
- Do not use surgical instruments to grasp the catheter.
- Do not use excessive force when inserting a catheter into a vessel.
- Ensure the catheter is thoroughly flushed and free of air prior to use.
- Avoid contact with liquids other than blood, saline, or contrast solution.
- Use in conjunction with fluoroscopic guidance and proper anticoagulation agents.

Flushing the catheter – Use the side port on the hub to flush the catheter. The catheter must be thoroughly flushed and free of air prior to use. When using the side port, consider sealing the proximal opening of the hemostatic valve with a thumb or forefinger.

5 Potential adverse events

- air embolism
- allergic reaction to contrast media
- arteriovenous fistula formation
- bleeding at the insertion site
- brachial plexus injury
- cardiac tamponade
- dislodgement
- dissection
- endocarditis
- heart block
- hematoma formation
- hemothorax
- infection
- irregular heart beat
- mediastinal widening
- perforation
- pneumothorax
- subclavian artery puncture
- thrombophlebitis
- thrombosis
- valve damage
- vascular occlusion
- vessel damage

6 Directions for use

6.1 Device preparation

1. Flush the C315 with saline.
2. Place the dilator in the C315.
3. Flush the dilator with saline.

6.2 Recommended procedure

1. Obtain venous access.
2. Advance the C315 and dilator together over a guidewire into the right atrium using fluoroscopy.
3. Remove the dilator.
4. Position the C315 catheter in the desired location.
5. Remove the guidewire.
6. Advance a lead or catheter through the C315 to the desired location.
7. When a lead has been positioned in the desired location, remove the C315 by slitting.
8. For instructions on removing the C315 by slitting, refer to the technical manual packaged with the slitter.

7 Detailed device description

7.1 Specifications (nominal)

Parameter	Models C315
Guide-catheter dilator	Material: Polyethylene
	Inner diameter: 0.96 mm (2.6 Fr, 0.038 in) minimum
	Outer diameter: 1.75 mm (5.2 Fr, 0.068 in) maximum
Guide catheter	Material: Polyether block amide,
	Inner diameter: 1.8 mm (5.4 Fr, 0.072 in) minimum
	Outer diameter: 2.4 mm (7.0 Fr, 0.092 in) maximum
Compatible Leads	1.4 mm (4.2 Fr, 0.055 in) or smaller
Compatible Inner Transvenous Devices	1.75 mm (5.2 Fr, 0.068 in) or smaller

Model number	H20	S4, S5, J	S10, H40	HIS
Guide catheter	20 cm	30 cm	40 cm	43 cm
usable length	(7.9 in)	(11.8 in)	(15.7 in)	(16.9 in)
Guide-catheter dilator	28.5 cm	38.5 cm	51.5 cm	51.5 cm
usable length	(11.2 in)	(15.2 in)	(20.3 in)	(20.3 in)

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

Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, Appendix Technical Manual for SelectSecure™ MRI SureScan™ model 3830 is true and unaltered copies of the original documents.

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This medical device is intended for use in clinical practice by qualified physicians and medical specialists.

MRI Conditions of Use

Patients and their implanted systems must be screened to meet the following requirements prior to undergoing an MRI:

- The patient has no implanted lead extenders, lead adaptors, or abandoned leads.
- The patient has no broken leads or leads with intermittent electrical contact, as confirmed by lead impedance history.
- The SureScan pacing system is implanted in the left or right pectoral region.
- The pace polarity parameters are set to Bipolar for programming MRI SureScan to On.
- The SureScan device is operating within the projected service life.
- For patients whose device will be programmed to an asynchronous pacing mode when the MRI SureScan mode is programmed to On, no diaphragmatic stimulation is present at a pacing output of 5.0 V and at a pulse width of 1.0 ms.

SureScan systems can be MRI-scanned with equipment that has the operating characteristics shown in Table 9. Continuous patient monitoring must be provided during the MRI scan as described in the Medtronic MRI Technical Manual.

Table 9: MRI Operating Characteristics

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Scanner type	Horizontal field, cylindrical bore, clinical system for hydrogen proton imaging
Scanner characteristics	- Static magnetic field of one of the following strengths: - 1.5 T - 3 T - Maximum spatial gradient of ≤ 20 T/m (2000 gauss/cm) - Gradient systems with maximum gradient slew rate performance per axis of ≤ 200 T/m/s
Scanner operation	1.5T — MRI radio frequency (RF) power — Normal Operating Mode. - The whole body averaged specific absorption rate (SAR) must be ≤ 2.0 W/kg. - The head SAR must be ≤ 3.2 W/kg. 3T — MRI RF power — First Level Controlled Operating Mode or Normal Operating Mode: B_{1+RMS}* must be ≤ 2.8 μT when the isocenter (center of the MRI bore) is inferior to the C7 vertebra. - Scans can be performed without restriction when the isocenter is at or superior to the C7 vertebra.

* B_{1+RMS} is a calculation of the amount of RF energy outputted by the scanner, whereas SAR is an estimate of the amount of RF energy absorbed by the patient.

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Pacing leads are insulated wires that deliver electrical impulses from the pulse generator to the cardiac tissue and sense cardiac depolarization. The SelectSecure leads are bipolar. Bipolar leads have two electrodes; one electrode is at the lead tip and one is located more proximally in a ring configuration, the cathode and anode, respectively. In the bipolar systems, a pulse generator creates an electrical pulse that creates a voltage field at both electrodes. The tip (cathode) electrode is designed to create a voltage field that depolarizes the local tissue, which initiates a cardiac contraction.

Medtronic

Medtronic Inc., USA

710 Medtronic Parkway N. E.

Minneapolis, MN 55432, USA

www.medtronic.com

**underlying the
work of the
medical products
and their features**

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Characteristics		SelectSecure™ MRI SureScan™ 3830
Type		Bipolar
Chamber		Atrium/Ventricle
Fixation		Nonretractable screw-in
Length		59 cm, 69 cm, 74 cm
Connector		IS-1 BI
Electrode surface area	Helix:	3.6 mm ²
	Ring:	16.9 mm ²
Tip to ring spacing		9.0 mm
Lead body diameter		1.4 mm (4.1 Fr)
Catheter Introduction size		3.0 mm (9.0 Fr)
Helix length (exposed)		1.8 mm
Resistance	Unipolar	25 ± 5 Ω (59 cm) 29 ± 6 Ω (69 cm) 31 ± 6 Ω (74 cm)
	Bipolar	84 ± 19 Ω (59 cm) 99 ± 22 Ω (69 cm) 107 ± 16 Ω (74 cm)
Steroid		Beclomethasone dipropionate
Amount of steroid (target dose)		100.0 µg max
Mass, g		59 cm – 2.28±5% 69 cm – 2.60±5% 74 cm – 2.77±5%
Anchoring sleeve length, mm		21.85 (-0.38/+0.25)
Compatibility with MRI		Yes
The mean bipolar voltage threshold at a pulse width of 0.5 ms, V		< 2.5
The mean bipolar pacing impedance at a voltage of 5.0 or 2.5 V and a pulse width of 0.5 ms, Ohm		350-2000
DC Resistance, Ohm		59 cm: Inner Circuit (tip to helix) 25 ± 5 Ω (Ohm) Outer Circuit (tip ring to distal ring) 59 ± 14 Ω (Ohm) 69 cm: Inner Circuit (tip to helix) 29 ± 6 Ω (Ohm)

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	Outer Circuit (tip ring to distal ring) $70 \pm 16 \Omega$ (Ohm) 74 cm: Inner Circuit (tip to helix) $31 \pm 6 \Omega$ (Ohm) Outer Circuit (tip ring to distal ring) $76 \pm 10 \Omega$ (Ohm)
Composite Tensile Strength, lbf (N)	≥ 1.0 (4.45)
IS-1 Pin to Connector Sleeve Tensile Strength, lbf (N)	≥ 2.5 (11.12)
IS-1 Connector Flex (dry)	Flex to 82.000 cycles at approximately 2 Hz with no intermittency or conductor fractures.

All technical specifications are subject to tolerances of $\pm 5\%$, unless otherwise specified.

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13. Viva XT CRT-D, model DTBA2D1, Viva S CRT-D, model DTBB2D1, Viva XT CRT-D, model DTBA2D4, Viva S CRT-D, model DTBB2D4, Viva Quad XT CRT-D, model DTBA2QQ, Viva Quad S CRT-D, model DTBB2QQ (RC # P3H 2017/6515 of 27.11.2017);

14. Implantable cardioverter defibrillator CLARIA MRI Quad CRT-D SureScan, model DTMA2QQ, Implantable cardioverter defibrillator CLARIA MRI CRT-D SureScan, model DTMA2D4, Implantable cardioverter defibrillator AMPLIA MRI Quad CRT-D, Implantable cardioverter defibrillator AMPLIA MRI CRT-D SureScan, model DTMB2D4, Implantable cardioverter defibrillator COMPIA MRI Quad CRT-D SureScan, model DTMC2QQ, Implantable cardioverter defibrillator COMPIA MRI CRT-D SureScan, model DTMC2D4, Implantable cardioverter defibrillator AMPLIA MRI Quad CRT-D SureScan, model DTMB2Q1, Implantable cardioverter defibrillator AMPLIA MRI CRT-D SureScan, model DTMB2D1, Implantable cardioverter defibrillator CLARIA MRI Quad CRT-D SureScan, model DTMA2Q1, Implantable cardioverter defibrillator CLARIA MRI CRT-D SureScan, model DTMA2D1, Implantable cardioverter defibrillator COMPIA MRI CRT-D SureScan, model DTMC2D1. (RC # P3H 2021/14463 of 28.05.2021).

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Disclaimer of warranty and general warning

Disclaimer of warranty

Although the Medtronic® Accessory, hereafter referred to as “Product” has been carefully designed, manufactured and tested prior to sale, the Product may fail to perform its intended function satisfactorily for a variety of reasons. The warnings contained in the Product labeling provide more detailed information and are considered an integral part of this DISCLAIMER OF WARRANTY. Medtronic, therefore, disclaims all warranties, both express and implied, with respect to the Product. Medtronic shall not be liable for any incidental or consequential damages caused by any use, defect or failure of the Product, whether the claim is based on warranty, contract, tort or otherwise.

The exclusions and limitations set out above are not intended to, and should not be construed so as to, contravene mandatory provisions of applicable law. If any part or term of this DISCLAIMER OF WARRANTY is held by any court of competent jurisdiction to be illegal, unenforceable or in conflict with applicable law, the validity of the remaining portion of the DISCLAIMER OF WARRANTY shall not be affected, and all rights and obligations shall be construed and enforced as if this DISCLAIMER OF WARRANTY did not contain the particular part or term held to be invalid.

General warning

The Products in this package are designed for use with other Medtronic products. Despite the exercise of all due care in design, component selection, manufacture, and testing prior to sale, the Products may be easily damaged by improper handling or use because of their unavoidably fragile character, which is dictated by the unusual requirements of their application. Consequently, no representation or warranty is made that failure or cessation of function of the Products will not occur.

The warranty period of storage of the device is not more than 24 months (2 years) from the date of packaging of the sterile device.

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
Declaration Letter

Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby or affirm that the attached documentation, Appendix instructions for use for delivery catheter C315H20, C315H40, C315HIS, C315J, C315S4, C315S5, C315S10 is true and unaltered copies of the original documents.

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Maximum pressure without leakage, kPa	300-320
Mass, g	2.9±5%
C315S5	
Usable Length, cm	30 (+1.0 / -1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)
Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)
Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	2.89±5%
C315S10	
Usable Length, cm	40 (+1.0 / -1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)
Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)
Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	3.0±5%
C315J	
Usable Length, cm	30 (+1.0 / -1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)

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Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)
Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	3.0±5%
C315H20	
Usable Length, cm	20 (+1.0 / -1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)
Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)
Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	2.7±5%
C315HIS	
Usable Length, cm	43 (+1.0/-1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)
Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)

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Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	3.06±5%
C315H40	
Usable Length, cm	40 (+1.0 / -1.6)
Inner diameter	1.8 mm (min.)
Outer diameter	2.4 mm (max)
Compatible Leads	1.4 mm or smaller
Compatible Inner Transvenous Devices	1.75 mm or smaller
Hub/Shaft Tensile Strength, N	≥15
Shaft Tensile Strength, N	≥15
Soft Tip Tensile Strength, N	≥5
Curve Stiffness, gf (N)	≤148 (1.45)
Straight Shaft Distal Segment Compression, gf (N)	≤296 (2.90)
Maximum pressure without leakage, kPa	300-320
Mass, g	3.08±5%

All technical specifications are subject to tolerances of ±5%, unless otherwise specified.

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© 2016 Medtronic
The following are trademarks of Medtronic: Attain, Medtronic, SelectSite