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Asclepion
Laser Technologies

EG-Konformitätserklärung *EC Declaration of Conformity*

Asclepion Laser Technologies GmbH

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07747 Jena

Germany

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Essential requirements of Directive 93/42/EEC, Annex I, with all relevant harmonized standards

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with the requirements of the Medical Devices 93/42/EEC Annex II (registration no. 1901839-006-000 of the notified body "LGA – InterCert GmbH", Tillystr. 2, 90431 Nürnberg, no. 1275). ¹⁾ Any modification to the product, not authorized by Asclepion Laser Technologies GmbH, will invalidate this declaration.

Jena, May 21th 2013

Asclepion Laser Technologies GmbH

I.V.

Pietro Amatulli
Technical Director

FB 7331

Revision: A

Asclepion
Laser Technologies

For:

Development

Project number

Name of the p

Item number

QuadroStarPRO

List of devices and accessories

DEVICES	PART NUMBER
QuadroStarPRO ^{GREEN}	1490
QuadroStarPRO ^{GREEN} with scanner and skin cooling	1491
QuadroStarPRO ^{YELLOW}	1492
QuadroStarPRO ^{YELLOW} with scanner and skin cooling	1493
QuadroStarPRO ⁹⁴⁰	1494
QuadroStarPRO ⁹⁸⁰	1495
APPLICATION PARTS for QuadroStarPRO ^{Green} & QuadroStarPRO ^{Yellow}	PART NUMBER
Scanner handpiece with skin cooling (Standard for 1491 and 1493)	4048
Set focusing handpiece with 1 mm spot (Standard for 1490 to 1493)	4042
optional available spot size optics 0.5mm	4043
optional available spot size optics 1.5mm	4045
optional available spot size optics 2.8mm	4046
APPLICATION PARTS for QuadroStarPRO ⁹⁴⁰ & QuadroStarPRO ⁹⁸⁰	PART NUMBER
Set focusing handpiece with 3,5 mm spot	4057
LipoSet	4228
Endo Set NIR Q600	2836
Endo Fiber NIR Q 600	2835
Bare Fiber NIR Q 300L	2864
Bare Fiber NIR Q 600L	2865
APPLICATION PARTS for QuadroStarPRO	PART NUMBER
Set Bare Fiber handpiece	4213
Bare Fiber handpiece	4212
Bare Fiber R Q 200	2963
Bare Fiber R Q 400	2969
Bare Fiber R Q 600	2970

Prepared:

Vs.	Change
02	Ergänzu 808nm g
01	First ex

Check and r

Checked by

Released by

Grund.

FB ers

17.07.

Datum

FB 7331	List of safety relevant materials and components	Asclepion Laser Technologies
Revision: A	Identification No.: MMo/13/10/10-01 Identification code/year/month/day-number	Page 1 of 4 FB 7331 LoCC QuPRO V02_MMo131010-01.doc

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
Check and release:

	Division	Date	Name	Sign
Checked by		Nov 17, 2013	I. Streit	
Released by		18.11.13	Thieme	

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FB 7331	List of safety relevant materials and components	 Asclepion Laser Technologies
Revision: A	Identification No.: MMo/13/10/10-01 Identification code/year/month/day-number	Page 2 of 4 FB 7331 LoCC QuPRO V02_MMo131010-01.doc

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
1.	Electric supply					
1.1.	Power cord (Netzkabel) 550704001	Feller	100-240 VAC max. 6,35A	250 V 10 A	VDE, CE EN 60320-1	outside
1.2.	KMF Power Entry Module with Line Filter (Gerätestecker-Kombi-element inkl. Filter) KMF 1.1193.1 u. 4301.1403 550901086	Schurter	100-240 VAC max. 6,35A	125 VAC / 60 Hz 250 VAC / 50 Hz 10 A	VDE 126123 CE EN 60320-1 EN 60939-2, CCC, CSA, UL (US/CDN) E72928, RU	On the back
1.2.1.	Line Fuse 550901090	Schurter	T 6,3A H 250V ~	T 6,3A H 250V ~	E oder H	Part of pos. 1.2
1.3.	Distribution line (Klemmblock) 149001017					
1.3.1.	Terminal Block 788-304	WAGO	230 V 3A/1,5kV	250 V 3A/ 4kV	KEMA, UL E45172, CSA,VDE	Part of pos. 1.3
1.3.2.	Relais RT314024 788-154	Tyco/Schrack	230 V 8A/4kV	250 V 16A/ 5kV	VDE 6106, UL E45172, CSA	Part of pos. 1.3
1.4.	Med. Netzteil 12/24V RL0402DE-1224 550399109	Ultra Level Tech	100-240V 6.35- 3.25A	100-240V 6.35- 3.25 A	UL E257331, IEC60601-1 (DE02015163	inside
1.5.	FN2010b-6-06 551201030	Schaffner	100-240 VAC max. 6,35A	125 VAC / 60 Hz 250 VAC / 50 Hz 10 A	VDE 126123 CE EN 60320-1 EN 60939-2, CCC, CSA, UL (US/CDN) E72928, RU	On the back
2.	Laser					
2.1.	Lasermodul Genesis MX 577-5000 MTM 551101143	Coherent	5W, 577nm	5W, 577nm	Verified by ALT	inside
2.2.	Lasermodul Genesis MX 532-8000 MTM 551101141	Coherent	8W,532nm	8W,532nm	Verified by ALT	inside
2.3.	940nm QD-C6830-FH(4) 551101145	Quantel	30W, 940nm	30W, 940nm	Verified by ALT	inside
2.4.	980nm QD-C6830-FH(4) 551101146	Quantel	30W, 980nm	30W, 980nm	Verified by ALT	inside
2.5.	Module M1F-SS5.2 551101152	DILAS	30W, 980nm	30W, 980nm	Verified by ALT	inside

Revision: A	Identification No.: MMo/13/10/10-01 Identification code/year/month/day-number	Page 3 of 4 FB 7331 LoCC QuPRO V02_MMo131010-01.doc
-------------	---	---

place of assembly

outside

On the back

Part of pos. 1.2

Part of pos. 1.3

Part of pos. 1.3

inside

On the back

inside

inside

inside

inside

inside

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
3.	Switch on circuit					
3.1.	PCB SOC 7601 06 032					Inside
3.1.1.	MP05-xx ME05-xx	Cotek	Input 230 VAC, 0.5A Output DC 24V, 0,21A	Input 85-264 VAC, 0.13A Output DC 24V, 0,21A	UL60601-1 E324816 CUL 60601-1 CB TÜV IEC 60601-1	Part of pos. 3.1 PCB SOC
3.1.2.	Fuseholder 0031.8201, 0853.0551	Schurter	230 V 6 A	250 V 10 A	IEC60127, VDE, CSA, UL E39328	Part of pos. 3.1 PCB SOC
3.2.	Schaltelement Oktron: Leucht- Tastkontaktgeber BZLOO 1.005.5173 550101073 Key Switch BZOO/1.000.7706 550101040	Schlegel	12 V DC 1 A	250 VAC 6 A	EN 61058 IEC60947-5-1 TÜV RU, CSA, UL E72889	inside
3.3.	Foot switch 660806008 MKF 1S D 1S-MED SK12 Alternativ: 660806009 MKF 1SD1S-MED GP12	Steute	5 VAC < 20 mA	max. 25 VAC max. 1 A	EN 60601-1 IEC 601 testet by RW TÜV IPX8, CSA, UL E178188	outside
4.	Controll system					
4.1.	Control Unit 760106066	ALT/ACD	+12 V < 2 A	+8-40 V 4 A	UL, Verified by ALT	inside
4.1.1.	Circuit board 760106066	ALT/ACD			UL E178188	Part of pos. 4.1
4.1.2.	DC/DC Converter RKZ 1205S	ReCOM	I:12Vdc/<0.2A O:5Vdc/<60m A Iso: 1.5 kV	In: 12Vdc/490mA Out: 5Vdc/1.2A Iso: 3 kV	EN 60601-1 60950-1	Part of pos. 4.1
4.1.3.	Digital Isolator ADuM2400ARWZ	Analog Devices	5Vdc/≤60 mA 4 kV	7Vdc/≤80 mA 4Kv	UL E214100, CSA60950 CSA61010 CSA60601 IEC60601 VDE 0884 TÜV6090 TÜV61010	Part of pos. 4.1
4.1.4.	Relais ALD 112	Panasonic	+5 V/< 40 mA 4 kV	277VAC/3A 30VDC/3A 4 kV	UL E43028, CSA, VDE EN 61810-1	Part of pos. 4.1
4.2.	Service-USB-Platine connection board 760106061	ALT/ACD	Input 12V, 0,5A	nput 12V, 0,5A	Verified by ALT	inside

FB 7331	List of safety relevant materials and components	
Revision: A	Identification No.: MMo/13/10/10-01 Identification code/year/month/day-number	Page 4 of 4 FB 7331 LoCC QuPRO V02_MMo131010-01.doc

no	type and item number	manufacturer	operating ratings	nominal ratings	approval	place of assembly
4.2.1.	DC/DC Converter TEN 5-1211	TRACO Power	I: 9-18V O: 5V/1A 1.5 kV	I: 12V O: 5V/1A 3 kV	CSA/UL 60950-1	Part of pos. 4.2
4.2.2.	Transformer 749014013	WE (Würth)	3.3 V 1.5 kV	3.3 V 4 kV	Des. IEC 60601-1 60950-1	Part of pos. 4.2
4.2.3.	ADUM 2201	Analog Devices	5 V 1.5 kV	5 V 5 kV	UL(E214100) IEC/CSA 60601-1 60950-1 61010-1 (Cnr:1646433) VDE0884-10 (Reg.nr:133860)	Part of pos. 4.2
4.3.	Glas touch unit 551102038	DOLD	+24 V 1 A	+24 V 1 A	Verified by ALT	outside
5.	Cooling System					
5.1.	peristaltic pump 66 0302 019	KNF Neuberger	12 V 12 W	12 V 12 W	CE	inside
5.2.	pressure switch 661202014	Suco	5V <0,1A	42V 4A	Verified by ALT	inside
5.3.	FAN AD0612XB- A72GL(P) 660301036	ADDA	12 V 0,375 A	12 V 0,375 A	CE	inside
5.4.	FAN 8212NG 660301032	ADDA	12 V 0,16 A	12 V 0,16 A	CE	inside
6.	Wiring					
6.1.	PVC- insulation tube 550806039, 550806043	Alpha Wire Coop.	2,5kV	15kV 4kV	UL /CSA E38155 Tested by ALT	inside
6.2.	wire 550702067	Lapp/ Helu	100-240 V	300/500V	HAR(H05) UL/ CSA	inside
6.3.	wire 550701117, 550701145	Lapp	100-240 V	300/500V	HAR (H05-V) UL/ CSA	inside
6.4.	Wire 550701141, 550701131	Lapp, Medikabel	100-240 V	300V	UL/ CSA	inside
7.	Enclosure					
7.1.	Plastic device cover Handpiece cover Biresin RG53 FR 149006002, 149006001	SiKa Deutschland	Device inside < 50 °C	110°C	UL 94, V-0 at 3 mm RU E204929	outside

User I
Ql

QuadroStarPRO

The laser device **QuadroStarPRO** emits high levels of emission in the visible or in the infrared range that may cause serious damage to the eye if the device is used improperly. Moreover, further risks are possible.

Therefore, this User Manual must be studied carefully before using the device and all warnings and instructions are to be observed during its operation.

It is of particular importance that all persons in the room in which the laser is operated put on the appropriate eye protection during the active laser operation.

Asclepion Laser Technologies GmbH

Brüsseler Str. 10

D-07747 JENA

Germany

Tel.: +49 (0) 3641 / 7700 - 401

Fax: +49 (0) 3641 / 7700 - 402

E-Mail: service@asclepion.com

Table of contents

1	Introduction	1
1.1	General information	1
1.2	Symbols in the user manual	1
1.3	Intended use	2
1.3.1	Indications	2
1.3.1.1	<i>QuadroStarPRO^{YELLOW}</i> (wavelength 577nm) and <i>QuadroStarPRO^{GREEN}</i> (532nm)	2
1.3.1.2	<i>QuadroStarPRO⁹⁴⁰</i> (940nm) and <i>QuadroStarPRO⁹⁸⁰</i> (980nm)	3
1.3.2	Contraindications	3
1.3.2.1	<i>QuadroStarPRO^{YELLOW}</i> (577nm) and <i>QuadroStarPRO^{GREEN}</i> (532nm)	3
1.3.2.2	<i>QuadroStarPRO⁹⁴⁰</i> (940nm) and <i>QuadroStarPRO⁹⁸⁰</i> (980nm)	4
1.3.3	Restriction of use	4
1.4	Declaration of Conformity as Medical Product	4
1.5	Copyright	4
2	Safety	5
2.1	General information	5
2.2	Electrical hazards	5
2.3	Biological hazards	6
2.4	Radio interference	6
2.5	Laser	7
2.5.1	Optical hazards caused by laser emission	7
2.5.2	Laser induced fire hazard	9
3	Technical Data	10
3.1	Product versions	10
3.2	Specifications	10
3.3	Technical description	12
3.4	Safety units	13
4	Start-up	16
4.1	Scope of delivery	16
4.2	Installation	17
4.2.1	Device design	17
4.2.2	Labels at the device	19

4.2.3	General information about installation	22	6.3.10
4.2.4	Door interlock connection	22	6.3.11
4.2.5	Foot switch	23	6.3.12
4.2.6	Optical fibres	24	6.4
4.2.7	Focusing handpiece (for QuadroStarPRO^{YELLOW} and QuadroStarPRO^{GREEN})	24	6.4.1
4.2.8	Scanner with integrated skin cooling (optional for QuadroStarPRO ^{YELLOW / GREEN}).....	25	6.4.2
4.2.9	Mains connection	25	6.4.3
4.3	Switching On	26	6.4.4
5	Operation of the device	27	7
5.1	General information on operation.....	27	
5.2	TOOLS menu.....	28	8
5.3	Adjusting the modes.....	29	8.1
5.3.1	BASIC mode.....	29	8.2
5.3.2	EXPERT mode.....	30	8.3
5.3.3	BURST mode.....	30	
5.3.4	SCAN mode.....	30	8.4
5.3.5	CW mode.....	31	
5.3.6	VARIX mode.....	31	9
5.4	Starting the treatment.....	32	9.1
5.4.1	Replacement of the handpiece, scanner or fibre.....	33	9.2
5.4.2	Switching off the device.....	33	
6	Application	34	10
6.1	General instructions.....	34	11
6.2	General treatment information.....	35	11.1
6.3	Skin treatment with the QuadroStarPRO ^{YELLOW} and the QuadroStarPRO ^{GREEN}	35	11.2
6.3.1	Interaction of the laser beam and the tissue.....	36	12
6.3.2	Classification of skin types.....	38	
6.3.3	Application parameters.....	39	13
6.3.4	Application techniques.....	40	13.1
6.3.5	Selection of patients.....	41	13.2
6.3.6	Contraindications.....	42	13.3
6.3.7	Possible side effects.....	42	
6.3.8	Preparing laser treatment.....	43	14
6.3.9	Test treatment.....	43	

.....22	6.3.10	Treatment procedure	45
.....22	6.3.11	Post treatment	45
.....23	6.3.12	Follow up treatment	46
.....24	6.4	Treatment with the infrared version of the <i>QuadroStarPRO</i>	47
GREEN)24	6.4.1	Interaction of the NIR laser beam and the tissue	48
LOW / GREEN)25	6.4.2	Application parameters	50
.....25	6.4.3	Application techniques	50
.....26	6.4.4	Selection of the patients, contraindications, side effects	54
.....27	7	Accessories	59
.....27			
.....28	8	Cleaning, disinfection and sterilisation	60
.....29	8.1	General rules	60
.....29	8.2	General advice for processing of the application parts	61
.....30	8.3	Special recommendations for cleaning and disinfecting the handpieces and scanners	64
.....30	8.4	Special recommendations for cleaning, disinfecting and sterilising the bare fibre and other special application sets	66
.....31			
.....31	9	Fault messages and fault removal	67
.....32	9.1	General information	67
.....33	9.2	Fault messages and fault removal	68
.....33			
.....34	10	Customer service	70
.....34			
.....35	11	Regular maintenance, safety checks and calibration	71
.....35	11.1	Routine maintenance	71
.....35	11.2	Regular safety check	72
.....36			
.....38	12	Disposal	73
.....39			
.....40	13	Manufacturer's declaration of EMC compliance	74
.....41	13.1	Emitted electromagnetic interference	75
.....42	13.2	Electromagnetic Immunity	76
.....42	13.3	Recommended Safety Distances	78
.....43			
.....43	14	EC declaration of conformity as medical device - Copy	79

15	Appendix: Cosmetic treatment with the infrared QuadroStarPRO.....	80
15.1	Intended use.....	80
15.2	Changes in other chapters	80

1 Int

1.1 Ge

**Kr
m:
sa**

**1.2 Syr
In
in**

1 Introduction

1.1 General information

- This User Manual contains the information required for the intended use of the device.
- It is an integral part of the complete product according to the applicable national, European and international directives and laws on product liability.
- The most important objective is to protect persons against risk situations and the device against damages caused by improper use.
- Regular maintenance according to the instructions given in this manual increases the precision and operability of the device during its entire service life.
- In case of malfunctions or errors that go beyond the scope of intervention of the User contact the Customer Service of Asclepion (or authorized distributor).

Knowledge of this User Manual is required for the operation of the devices. Therefore, please make yourself familiar with the contents and pay special attention to information concerning the safe operation of the device.

1.2 Symbols in the user manual

In this manual, the following symbols are used to draw your attention to dangers or to give instructions on operating the device:

DANGER

This symbol indicates a direct danger for the life and health of persons. Non-compliance with these instructions causes serious danger to health up to life-threatening injuries or even death.

This symbol indicates a possible danger for the life and health of persons. Non-compliance with these instructions can cause serious danger to health up to life-threatening injuries.

CAUTION

This symbol indicates a possible dangerous situation. Non-compliance with these instructions can cause injuries.

Important information

This symbol indicates important information about the proper use of the device. Non-compliance with these instructions can cause damage to property or malfunctions of the device or in its surroundings.

Tip

This symbol indicates practical tips and particularly useful information that enables the user to use all the features of the device in an optimal manner.

1.3 Intended use

The incorrect using of this device can cause serious danger to health. Therefore, all warnings, instructions and information of this User Manual are to be observed under all circumstances! The non-compliance is considered as unintended use of the device!

Any modification of this device is not allowed!

The *QuadroStarPRO* is a laser (Class 4 laser product) that has relatively long impulses up to a continuing emission at all varieties of wavelengths (see chapter *Technical Data*). Therefore the device can cause a selective thermal impact or damage to the tissue and hence can be used for various indications. The versions in the visible (green and yellow) and in the infrared range are nearly identical in their effect. For this reason the intended use and indications are given together for the specific wavelength.

1.3.1 Indications

1.3.1.1 *QuadroStarPRO*^{YELLOW} (wavelength 577nm) and *QuadroStarPRO*^{GREEN} (532nm)

The wavelengths in the green and yellow range are very well absorbed by melanin and haemoglobin. Hence the laser is ideal for the treatment of superficial vascular lesions. The gentle removal of benign pigmented lesions and other skin changes by means of coagulation and vaporisation is also gently possible with these wavelengths.

In detail the following skin lesions are treated for many years now and are by now standard procedures^{[1], [2]}.

- Teleangiectasis (rosacea, couperosis, rubeosis faciei, Erythrosis interfollicularis, Rendu-Osler-Weber disease) also associated with other diseases (lip angioma, venous angioma, NBCCS)
- Small granuloma telangiectaticum (pyogenicum)
- Spider naevi
- Cherry haemangioma (Campbell de Morgan spots)
- Angiofibroma (Brookes disease, Bourneville-Pringle's disease)
- Extensive, plane haemangioma
- Extensive, plane naevi flammei
- Spider veins (red, small and superficial vessels)
- Xanthelasma
- Syringioma
- Warts
- Poikiloderma civatte
- Skin changes due to Kaposi's sarcoma
- Adenoma Sebaceum
- Chondrodermatitis nodularis heliis
- Benign epidermal pigmented lesions
- Striae atrophicae distensae (improvement, decrease of the reddening but NO removal)

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[2] Vascular lasers and IPLs: Guidelines for care from the European Society for Laser Dermatology (ESLD)
M.Adamic et al., Journal of Cosmetic and Laser Therapy, 2007;9:113-124

1.3.1.2 On

Wa
rea
for
Dis
the
dur
De
the
do
In
giv
For
for

1.3.2 Co

1.3.2.1 On

1.3.1.2 QuadroStarPRO⁹⁴⁰ (940nm) and QuadroStarPRO⁹⁸⁰ (980nm)

Wavelengths near the infrared range (NIR) are absorbed by haemoglobin and water, yet they reach deeper into the tissue than yellow or green laser beams, due to the relatively low absorption and scattering.

Diode lasers in the NIR have a variety of indications in surgery and the endoscopy. Yet, due to the considerably lower selectivity in comparison with the green and yellow laser beams (absorption of haemoglobin and water), these devices are usually not suitable for dermatological or aesthetic indications, except for very general ones like warts or discoloured and/or thickened nails due to Onychomycosis (see attachment).

Depending on the fluence and power density the soft tissue can be coagulated or vaporised. In the fiber contact method the tissue can also be cut, as this is a special case of vaporisation. The diode lasers are also very suitable for endovascular coalition of enlarged veins (varicose veins).

In the dental sector the diode lasers are used in the soft tissue surgery, as with the application given above and for dental bleaching.

For more details and Information about possible indications and applications see chapter *Application / Treatment with the infrared version of the QuadroStarPRO*.

1.3.2 Contraindications

CAUTION

The intensive laser light in the visible or infrared range can cause intentional and unintentional thermal damage of skin and tissue structures. The contraindications listed below are based on clinical experience in operating laser devices for said purpose during many years. They do not claim to completeness or unlimited validity.

The definite decision on the patients' suitability for the individual treatment in question is the sole responsibility of the attending physician.

1.3.2.1 QuadroStarPRO^{YELLOW} (577nm) and QuadroStarPRO^{GREEN} (532nm)

The intensive laser light in the infrared range can cause thermal damage to skin structures. Therefore it is **NOT PERMITTED to use laser treatment** for pregnant women or children because clinical studies do not exist for these clients.

- Tanned skin (postpone treatment)
- Sensitivity to light, usage of preparations that increases sensitivity to light
- Dark skin type (IV to VI)
- Cancerous and precancerous lesions and naevus cell naevi in the treatment area
- Kasabach-Merrit-Syndrome (haemangioma with coagulopathy)

The (unintentional) **laser treatment** of precancerous lesions, malignant skin changes or of the Kasabach-Merrit-Syndrome poses a **significant risk for the health of the patients**.

It is an absolute prerequisite that in case of the slightest doubt about the benignancy of pigmentation or any form of haemangioma, the diagnosis has to be confirmed or treatment has to be cancelled.

- Herpes simplex in the treatment area (prophylaxis or postpone treatment 10 days)
- Tendency to excessive scar formation (therapy of keloids in combination with other laser devices possible)

1.3.2.2 QuadroStarPRO⁹⁴⁰ (940nm) and QuadroStarPRO⁹⁸⁰ (980nm)

As these devices can be used for various applications in different specialist discipline, the assessment of contraindications is the solely responsibility of the attending physician.

At the same time it has to be borne in mind, that especially the use of optical fibres in a non contact method and during a long emission can cause a coagulation of several millimetres in the tissue, which can be fatal for the patient.

1.3.3 Restriction of use

The laser emission of the device can cause thermal changes of tissue structures.

Treatment may not be performed by anyone other than a practising physician (referred to as user for simplification hereafter)

Active medical products may only be applied by users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience. For skin treatment in Germany the user is advised to gain a qualification verification

It is expressly pointed out that the use of the device is only permitted for Users that have completed successfully a laser safety course including application training and have been instructed to operate the device. This has to be confirmed with a signature in the medical device book.

Asclepion Laser Technologies GmbH also offers these courses for laser safety in medical application.

1.4 Declaration of Conformity as Medical Product

The QuadroStarPRO has been designed to comply with the EC Directive 93/42/EEC regarding medical products and its national implementation in the form of the German Medical Product Act (MPG). This conformity declaration refers to the medical intended use given in chapter 1.3 solely. It is not to apply to the cosmetic application given in the appendix.

The CE marking CE 1275 refers to the device for medical intended use solely.

Device class under Directive 93/42/EEC:

II b

UMDNS-No.:

16-948 Laser, therapeutically

1.5 Copyright

The specifications are subject to change due to further technical development. Contact the Technical Service of Asclepion Laser Technologies GmbH or your distributor of Asclepion Laser Technologies GmbH to get latest information.

© Unless expressly authorized, the forwarding and duplication of this document, and the utilization and disclosure of its contents are not permitted. Violations will entail an obligation to pay compensation.

All rights reserved in the event of granting of patents or registration of a utility model.

2 Safe

2.1 Gen

The 93/42

In ac
devic
invas

Like
be a
biolo

Obst
devic
the r

2.2 Ele

2 Safety

2.1 General information

The *QuadroStarPRO* complies with the requirements of the EC Directive for Medical Devices 93/42/EEG and the German Medical Product Act (MPG).

In accordance with this EC Directive and its Annex IX, the *QuadroStarPRO* qualifies as a class IIb device. (This is not to be confused with the laser class, as it defines the device as a surgically invasive product).

Like each device, the *QuadroStarPRO* may lead to potential risks when operated. The user should be aware of them before starting the system. Among these risks are optical, electrical and biological dangers as well as fire hazards.

Observe the national regulations applicable for the operation of this product. In particular, such devices may be operated only in accordance with the generally accepted rules of technology and the relevant regulations on labour safety and the prevention of accidents.

The laser device emits high levels of radiation in the visible or infrared range that may cause serious damage to the eye if the device is used improperly.

Absolutely observe the instructions given in this User Manual for equipment and procedures for operating, maintaining, checking or calibrating purposes.

The laser radiation of the device can cause, whenever low, the risk of fire. Under no circumstances the device may be used in an explosive atmosphere (classified as AP and APG under IEC 601-1). Take care to preclude the existence of vapours from solvents or flammable liquids (these have been possibly used for cleaning or disinfection) in the room before you start working with the laser.

Please note that this product is subject to scheduled technical safety tests that must be performed at annual intervals to ensure the safety for the staff and clients. If the tests are not performed regularly serious health problems can be caused.

Work for service, repair or modification may not be performed by anyone other than Asclepion Laser Technologies GmbH personnel or other bodies duly authorized to handle such work to ensure the safety for the staff and clients.

If unauthorized persons carry out service works, repairs or modifications at the device serious health problems can be caused.

2.2 Electrical hazards

The laser is operated with line voltage. There is the risk of an electrical shock! Do not open the device, if you are neither trained nor authorized to do so.

The device is only to be operated with reliably grounded and perfectly installed sockets with earthing contact.

The laser is operated with line voltage. There is the risk of an electrical shock!
The device is designed to prevent the penetration of liquids during normal use. If liquid penetrates into the inner system space or the handpiece in exceptional cases, you must turn power off and pull the main power plug. Please contact our Technical Service.

2.5 Laser

2.5.1 Optics

Lasers
class 4

2.3 Biological hazards

Laser fumes may contain viable tissue particles if tissue is vaporized. There is the risk of microorganisms being transferred and damage of the respiratory tracts due to the smoke itself. You are advised to provide a fume extraction device with filters for vaporization or ablation of tissue.

During treatments of the skin the contact surface of the handpiece may have contact with the patient's skin. When using fibres and other application parts for surgical and endoscopic purposes, contact with tissue and body fluids is usually unavoidable. There is a risk of microorganisms being transferred.

The used handpiece must be thoroughly cleaned and disinfected immediately after each application. The fibres or application parts for surgical and endoscopic procedures have to be disposed in case of disposable items or to be thoroughly cleaned, disinfected and sterilised (see chapter *Cleaning, Disinfection and Sterilisation*).

For fu
standa
beams
danger
applic
Please

It is hi

1.

2.

2.4 Radio interference

The *QuadroStarPRO* fulfils the requirements of standard EN 60601-1-2. The system is not influenced by electromagnetic radiation from other devices that meet the same standard, nor does it generate any electromagnetic radiation of a kind that exceeds the limiting values under EN 60601-1-2. A Manufacturer Declaration on Electromagnetic Compatibility issued per EN 60601-1-2 is attached to this User Manual.

Medical products are subject to special safety precautions regarding electromagnetic compatibility (EMC). The device emits low electromagnetic radiation, that may influence other devices not complying with the above given standard. The result may be unintended settings.

Please follow special instructions and advisory notes in the chapter *EMC Manufacturer Declaration for installation and operation of the QuadroStarPRO*.

The device may be influenced by devices that do not comply with the standard indicated above. The result may be unintended settings. Therefore, switch off cellular phones and similar equipment before putting the laser device into operation.

The use of accessories and cables not authorized by the manufacturer may result in an increased interference emission or a decreased immunity of the device with the above mentioned consequences.

Use only accessories and cables that have been authorized by Asclepion Laser Technologies for this device.

3

The device must not be stacked / arranged close to or in combination with other devices. If the device is operated close to or with other devices it must be checked if the device is working properly in this arrangement.

2.5 Laser

2.5.1 Optical hazards caused by laser emission

Lasers are classified in accordance with their potential for danger. The most dangerous class is class 4 (the least dangerous class 1).

The device is a laser of class 4. That means that the direct and even scattered laser radiation can cause serious damage to the eye if the device is operated improperly.

The protection measures described in the following must be strictly observed. It is of particular importance that all persons in the room in which the laser is operated put on the appropriate eye protection (see chapter *Technical Data*) during the active laser operation.

For fundamental rules on the handling of laser devices you are referred to the international standard EN 60825-1 and EN 60601-2-22 supplement 1: "Guidelines for the safe use of laser beams on humans". It is complemented by national regulations providing general protection from dangerous laser radiation. Their purpose is to protect operating personnel and clients during applications.

Please follow specific national regulations that may be valid for your location.

It is highly recommended to follow the rules specified below:

1. A laser protection adviser should be appointed.
2. The NOHD (Nominal Ocular Hazard Distance) of this laser system is very great so that the entire laser session room is to be regarded the laser area. **Warning lights** (preferably yellow) and a triangularly shaped yellow **laser warning label** must be provided at all access points (for US see below right). The lights have to glow when the laser is in operation.


Figure 1: Laser warning label (left: European; right: US style)

3. Personal eye protection

Laser radiation can cause serious damage to the eye.

Everyone present in the laser room during a treatment session must wear laser safety goggles, the clients an eye shield!

The laser safety goggles and the eye shield for clients must comply with the specifications defined in the technical data section!

4. Further precautionary measures

- Cover up windows and other openings of the treatment room in order to prevent unwanted emission of laser radiation.

- The walls in the treatment area must have a rough surface. Remove all mirrors or other reflecting objects from the treatment area.
- Restrict access to the treatment room to the client to be treated and to those persons who will assist and are trained in the handling of this laser device. All personnel working in the treatment room should have got a laser safety instruction.
- Make sure that personnel trained and assisting in a therapy session know how to shut down the laser in the case of an emergency.
- Point the active laser only onto the area that is intended for treatment.
- Refrain from looking directly into the outlet opening of the handpiece or that of the optical fibre, even if you wear laser safety goggles.
- Remove all metal objects such as watches, rings, bracelets or similar belongings from the working range and refrain from the use of reflecting instruments or other comparable materials where possible.
(Reflecting objects can interrupt the laser beam and guide it to another area than the intended one. Many surfaces, even the ones that seem to be mat, can strongly reflect the emission wavelength of the laser.)
- During surgical, endoscopic or similar procedures additional safety measures have to be taken. Special attention has to be paid to the fact that the end of the fibre always protrudes a few millimetres from the working channel of endoscopes, cannulas or others during laser operation. This is to ensure that, due to generated heat, no instruments are damaged or fire breaks out, which can damage the patient's health. The instruments have to be made of a material, that is suitable for laser treatments.
- For phases of non-operation, you should switch to STANDBY mode (in the STANDBY mode accidental release of the laser is prevented).
- Always remove the key from the key switch once power has been turned off and keep the key in a safe place.

5. Responsibility

The user who releases the laser emission is responsible for the laser safety. He/she has to ensure that all preventive measures have been taken **before** activating the laser. It is of particular importance that all persons in the room, in which the laser is operated, including the clients, put on the appropriate laser safety goggles or eye shield.

252 L

2.5.2 Laser induced fire hazard

Due to the laser radiation of the devices there is a risk of fire or explosion if the laser is used in the presence of flammable materials, solutions or gases or in an environment which is enriched with oxygen.

The Users should take the following precautions in order to prevent fires caused by laser radiation:

- Please expose only the proposed body areas to laser emission
- Keep clothes as far as possible away from the treatment zone.
- Some materials, which are saturated with oxygen (e.g. cotton), could be inflamed at the high temperatures, which develop under the intended use of the laser.
- Please keep a minimum of flammable materials in the treatment room only. If it is necessary for the treatment to use flammable materials, please moisten these.
- Please use non flammable substances for the preparation of the treatment and for the cleaning and disinfection of the handpieces directly before and during the treatment. If it is not possible to avoid usage of solvents such as alcohol or isopropyl alcohol directly before and during the treatment, the vapour should be removed by aeration before activation of the laser or at least wait for some time.
- Have always a small amount of water (e.g. a glas) at hand in the treatment room to be able to extinguish smouldering towels or small fires immediately.
- If any electrical device catches fire it must not be extinguished with water. Otherwise there is a high risk of an electrical shock. Therefore always have a small fire extinguisher suitable for electrical devices at hand.
- In case of an emergency immediately activate the Laser emergency STOP switch in order to immediately stop laser emission.
- Please use non-flammable substances for anaesthesia (incl. skin cooling), if at all.
- Avoid the use of oxidizing gases, e.g. nitrogen oxide (N₂O) and oxygen. Be especially careful when using oxygen. Oxygen increases both the intensity and the extent of fire.
- If you use endoscopes, trocars, cannulas or similar devices, always pay attention, that they are not inflammable and that the end of the fibre always protrudes the instrument at least 7 to 8 mm during laser emission!!
- In order to avoid a deflagration of endogen gases (e.g. Methane in the digestive tract), use locally-acting ventilation systems.

3 Technical Data

3.1 Product versions

QuadroStarPRO^{GREEN} Item no. 1490	Green laser with handpiece
QuadroStarPRO^{GREEN} Item no. 1491	Green laser with scanner and integrated skin cooling
QuadroStarPRO^{YELLOW} Item no. 1492	Yellow laser with handpiece
QuadroStarPRO^{YELLOW} Item no. 1493	Yellow laser with scanner and integrated skin cooling
QuadroStarPRO⁹⁴⁰ Item no. 1494	Laser with 940nm
QuadroStarPRO⁹⁸⁰ Item no. 1495	Laser with 980nm

3.2 Specifications

Device model	Table-top unit
Display	LC Display
Operator guidance	Operating element with capacitive buttons Optional with iPad, separate user manual
Cooling	Room air
Door interlock connector	5 VDC/ 20 mA
Laser warning lamp	Potential-free relay contact max. 24 V / 1 A (normally open contact)
Permissible ambient conditions	Temperature: 15 °C to 30 °C rel. Air humidity: max. 70 % (no condensation) Height: max. 3000m above sea level
Permissible transport and storage conditions	Upright position Temperature: 5°C to +45°C rel. Air humidity: 10 % to 95 % (no condensation)
Dimensions	Basic body: 420 x 380 x 190 mm (W x D x H)
Weight	ca. 12 kg
Power requirements	100 - 240VAC, 50/60Hz, max. 635VA @ 100V, max. 780VA @ 240V
Main fuse	fuse 5 x 20 T6,3A H 250V (6,3A, slow-blow, high breaking capacity)
Overvoltage category	II (IEC 60664-1)

Max. Peak of mains voltage	V MT peak = 2500V peak
Max. Nominal mains voltage	V MN r.m.s. = 300V r.m.s.
Classification under Directive 93/42/EEC	II b
Accessories	See chapter 7 <i>Accessories</i>
Operating mode	The device has been designed for continuous operation. In case a temperature is reached, which does not allow to work, the device switches automatically for a short time to cooling mode until a normal operating temperature is reached.

Laser specifications	
Laser type	Diode laser
Wavelength	QuadroStarPRO ^{GREEN} : 532nm
	QuadroStarPRO ^{YELLOW} : 577nm
	QuadroStarPRO ^{940/980} : 940nm or 980nm
Laser class	4
Required laser safety goggles (acc. to DIN EN 207: 2009)	QuadroStarPRO ^{GREEN} : 532nm D L4
	QuadroStarPRO ^{YELLOW} : 577nm D L4
	QuadroStarPRO ^{940/980} : 940nm/ 980nm D L4
Pulse length	1ms - 95sec and cw (depending on the mode)
Power	QuadroStarPRO ^{GREEN} : max. 8W
	QuadroStarPRO ^{YELLOW} : max. 5W
	QuadroStarPRO ^{940/980} : max. 25W
Pulse frequency	Max. 20 Hz (depending on the mode)
Spot sizes (Laser beam diameter)	handpiece: 1mm, optional 0.5mm, 1.5mm, 2.8mm Scanner: 1mm
Tolerance of output power	max. ±20 %
Laser beam mode	Multimode
Beam divergence at the hand-piece output (round angle, 1/e ²)	typ. 250 mrad @ Handpiece typ. 250 mrad @ Scanner typ. 294 mrad @ Fibre
Nominal ocular hazard distance (NOHD)	Handpiece: 4.5m Fibre: 5.2m
Degree of protection, application part	B

Laser type	Pilot laser (laser diode)
Wavelength	650 nm
Laser class	2
Laser power	<1 mW, adjustable
Operating mode	cw or flashing (ca. 1 Hz)

3.3 Technical description

Power supply unit

The integrated power supply units transform the line voltage in rectified low voltage for the electrical supply of the device.

Control

The output parameters and the safety systems of the device are monitored and controlled by a micro control system.

Operating modes

The device operates in both repetition pulsed mode and continuous wave mode.

The device is designed for continuous operation in the specified ambient conditions.

In case a temperature is reached, which does not allow to work, the device switches automatically for a short time to cooling mode until a normal operating temperature is reached.

Operation

The device is operated with capacitive control elements. The display works with a dot matrix. Depending on the operating mode different parameters can be modified, like spot size, fluence or power, pulse duration, pulse frequency and others (see chapter *Operation of the device*).

In addition to the selected parameters other data, such as the laser operating state (STANDBY or READY), and, if necessary, error plain-text messages are displayed.

Optionally the device may also be operated via an iPad. Still the button READY can only be activated on the touch screen of the *QuadroStarPRO*, in order to provide a last control over the modified parameters by iPad at the device.

Cooling

The device cooling is ensured by the indoor air. External cooling water is not required

To maintain the cooling effect, the device should only be operated in the specified ambient conditions. The ambient air should be dustless. At the openings provided for cooling the device, a minimum-distance has to be kept to air-tight-objects (e.g. the wall of the room). (See chapter *Installation*).

Laser

~~The laser of the QuadroStarPRO is a laser of the laser class 4. Depending~~ on the activated mode the laser can emit impulses in the millisecond range up to a continuously emission. The wavelength is in the green, yellow or near infrared range (see chapter *Technical Data*).

The laser beam emission is activated by pressing the foot switch.

Skin Cooling

When using the scanner an efficient skin cooling is indispensable. For this reason the scanner contains a metal probe including a sapphire window, which is cooled by an integrated peltier element. The temperature of the probe is monitored automatically.

3.4 Safety units

The device is equipped with several safety units that have been provided to avoid maloperation and the unintentional activation of the system. All the persons who operate the laser or assist during the treatment should make themselves familiar with these units.

Laser emergency STOP switch

The red Laser emergency STOP switch at the front of the device (see figure below) is for the purpose of switching off the system immediately in an emergency. It should only be activated in emergency situations, i.e. if the emission has to be interrupted immediately.

To turn off the system immediately press this button. After the elimination of the emergency the device can be switched on again by the key switch.

Note

The Laser emergency STOP switch should not be used to turn the system on and off in its normal condition.


Figure 2: Safety units at the front panel of the device

Key switch

The key switch is used to activate the system. Only authorized persons having access to the switch can start the system.

The laser device emits strong radiation in the visible and the infrared range. This radiation can cause serious eye injuries if the device is not properly used.

- The device must not be switched on or operated **by unauthorized personnel**. Remove the key always after having switched off the system and ensure that it is kept by authorized persons only.

Power switch

The power switch separates the device from the line voltage. If the device is not in use, the switch should be in the OUT (to press 0) position.


Figure 3: key switch and power switch at the rear panel of the device

Footswitch

The foot switch is an electric button that releases the laser emission if the device is in the READY operating state. The switching elements are twice because of safety reasons. Position the foot switch always close to the treatment area.


Figure 4: Footswitch – left: standard – right: optional

ared range. This
ed.
onnel.
sure that it is kept

in use, the switch

is in the READY
Position the foot

– July 30th, 2013

You may order an optional foot switch on which you can increase or decrease one parameter by the means of a left and right contact switch (INTENSITY or SPEED, can be selected in the TOOLS menu). Although the modification of the parameter is not possible when the footswitch is activated and the laser emits radiation.

The foot switch has always priority – which means, that the release of the switch interrupts the emission of radiation immediately, independent on the pulse duration pre-set at the display.

STANDBY mode

The STANDBY operating state prevents the unintentional or inadvertent activation of the laser. If the device is in the STANDBY operating state, laser emission cannot be activated. The user can only start the emission, if the READY button has been pressed.

The system switches to STANDBY in the following situations:

- after the first start of the device,
- if the system has been in the READY operating state for a longer period of time without the laser released,
- if the User presses the STANDBY key when the system is in the READY operating state.

Automatic shut-down

The device is equipped with an automatic shutdown system. If a specific problem arises, the system will automatically change into the following safe mode:

The laser beam shutter is closed, the laser discharge is interrupted and the footswitch is deactivated. The screen shows a fault message identifying the specific error (see chapter *Fault Messages and Fault Removal*).

Remote door interlock

The system is equipped with a remote interlock contact that may be connected to the entrances to the laser room. If the remote interlock contact is open (e.g. if the door is opened), the laser emission will be interrupted automatically and the device will change to the STANDBY mode.

Acoustic signal

The laser emission is indicated by an acoustic signal during the period of emission.

Warning labels

The device carries several warning labels (see chapter *Installation - Labels at the device*). These labels must always be clearly readable and replaced immediately when damaged.

Contact the Technical Service immediately in such cases.

The manufacturer will only guarantee the safety, reliability and performance of the device, if:

1. the device has been used in accordance with the instructions given in this manual,
2. the installation, assembly, extensions, changes, repair and maintenance works have been performed by persons being authorized and qualified for them, and
3. the electric system at the place of installation complies with the requirements placed by the EN / IEC and the local regulations.

4 Start-up

4.1 Scope of delivery

Important Information

Check together with the haulier whether all components that are listed below and required for the safe operation of the device are present and not damaged. Incompleteness and/or damages in transit are to be claimed immediately at the supplier.

Asclepion Laser Technologies does not accept any liability unless another agreement has been made expressly.

Please, keep the packing for the future safe consignment of the device.

CAUTION

Refrain from start-up action of any kind if the device or its accessories are found to have suffered mechanical damage. There is risk of physical injury if this advice is disregarded.

Please contact the Technical Service of Asclepion Laser Technologies GmbH in such cases.

The basic equipment comprises the following device parts:

- basic device
- power cable
- foot switch
- Application part (handpiece / scanner / fiber, dependent on device version)
- door interlock connector
- keys (two for key switch)
- laser safety goggle, eye protection for client
- user manual
- laser warning sign
- Set for refilling the water circuit (optional, only for scanners with skin cooling)

4.2 Install

4.2.1 Device

Error indicator

Laser emission
STOP switch

Handpiece on
scanner holder

Service

Fiber, SMA con

4.2 Installation

4.2.1 Device design


Figure 5: front side of the device



Figure 6: Connections at the rear side of the device

Connectors in detail:

	Power switch 1 = On (to press 1) 0 = Off (to press 0)
	Mains fuse Micro fuse 100-240 V~ 50/60 Hz 6,3A slow-blow, breaking capacity H, 5x 20
	Key switch  On  Off
	Footswitch connector
	2 ports: door interlock and door warning lamp
	Scanner, electrical connector
	Scanner, water connector (cooling circuit for skin cooling)
	Scanner, venting cooling circuit (skin cooling) Luer-Lok connector
	Optical fibre, SMA connector (optical)
	Optical fibre, electrical connector

422 Lb

4.2.2 Labels at the device

The laser device emits high levels of radiation in the visible and the infrared range that may cause serious damage to the eye if the device is used improperly.

Therefore, warning and information labels are attached at appropriate points of the device. All the staff must know these labels and their meaning to ensure the proper use of the device.


Figure 7: Labels at the front side






Figure 8: Labels at the rear side

Labels in detail:

	Laser emergency STOP switch This label indicates the Laser emergency STOP switch to immediately switch off the laser radiation.
	
	Laser warning labels EN 60 825-1 for <i>QuadroStarPRO^{GREEN}</i> (532nm) <i>QuadroStarPRO^{YELLOW}</i> (577nm) <i>QuadroStarPRO^{NIR}</i> (940, or 980nm) This label informs about the dangers, the maximum values of emission and the classification of the laser source.
	Attention: follow User manual: This label requires the user to strictly follow the requirements given in the User manual.
	Laser warning label EN 60 825-1: This label warns of the hazard of the device emitting laser radiation.
	Laser beam aperture : This label indicates the exit point of the laser beam out of the transmission system.
	Application part B

SERVICE

Manufacturer Model Serial number Date of manufacture: YYYY-MM CE symbol Protection degree application parts (B) Permissible operating temperature		 Asclepion Laser Technologies GmbH 07747 Jena, Brüsseler Straße 10, Germany Model / SN QuadroStarPRO 1490 / 00100 100-240 V~ 50/60 Hz 100 V~ 635 VA 240 V~ 780 VA  2012-11  	Type plate Indicating type-specific data of the device such as the serial number and the date of manufacture.
100-240V ~ 50/60Hz  T6.3A H 250V ~		Power supply This label indicates the mains requirements and the value of the mains fuse.	
		Footswitch connector	
		Door warning lamp connector	
		Door interlock connector	
		Scanner connector including skin cooling: Connectors for cooling water on the left and right and electrical plug in the middle	
		Venting opening; necessary when refilling water circuit of the skin cooling	
SERVICE    ① ②		Interface labels: SERVICE: Only for authorised service personnel ①: USB Interface ②: Emission of radio frequency/ WLAN	

Emergency STOP
User radiation.

the maximum
ion of the la-

following the re-

device emit-

the laser beam

July 30th, 2013

4.2.3 General information about installation

Important information

After the installation and any relocation of the device from a cold to a warm environment with a temperature difference of more than 5°C, allow the device to adjust to the room temperature in unpacked condition before being used at least in the following way:

min. 2 hours at a temperature difference of up to 10°C

min. 4 hours at a temperature difference of up to 15°C

min. 8 hours at a temperature difference of more than 20°C.

Afterwards, switch the device on without a handpiece connected and allow the device to warm up for at least 30 minutes.

Take care that the installation and operating environment of the device meets the following conditions: temperature: 15 to 30°C, relative humidity: max. 70 % and height not more than 3000 m above sea level.

Non-observance of these instructions may result in the destruction of the device.

For an efficient cooling of the device, please keep a minimum distance of at least 20 cm between its lateral and rear sides to the wall.

The laser device emits high levels of radiation in the visible and the infrared range that may cause serious damage to the eye if the device is used improperly.

Every entrance to the laser area (usually the laser treatment room) must be equipped with warning lamps that are lighting (or flashing) as long as the laser emits radiation to warn persons against entering the room.

4.2.4 Door interlock connection

At the rear panel of the device, a connector is provided for the connection of a door interlock contact to the device (see figure below).

As a port a round plug has to be used which can be screwed, type (x) V 40 according to IEC 60130-9, e.g. Lumberg WSV 40 (as delivered).

(x): dependent on type, e.g. WS for an angled plug as delivered.


Figure 9: Door Interlock and door warning lamp connectors

By connecting the two connections of the *door interlock contact* with an external door contact you can prevent laser radiation from being emitted when the door of the laser room is opened.

This safety system serves to shut off the treatment beam automatically when the door is being opened thus precluding any risk to the person entering.

The plug of the *door interlock contact* delivered is factory-fitted with a jumper. To connect a door interlock contact, remove this jumper and connect your door interlock contact to the same terminals. The external contact must be potential-free and designed for at least the current and voltage specified in the chapter *Technical Data*. If a door-interlock contact is not used, make sure to keep the plug with the factory-fitted jumper connected.

Important information

The contact plug may only be connected with the device being turned off and disconnected from the power supply. Non-observance of this instruction can cause damages of the device.

The contacts for the remote interlock and the foot switch must never be connected with line power because otherwise the system will be seriously damaged. Connect these two connectors only in the way described in this section.

Tip

If you have connected the *door interlock contact* with the entrance door, please check every day before starting a treatment that the laser cannot be switched to the **READY** with this door being open. In case of malfunction please inform the Technical Service.

An external voltage source for a low-voltage warning lamp can be switched via two contacts of the warning lamp connector. The internally installed switching relay is designed for a maximum voltage and a maximum current given in the chapter *Technical Data*. Absolutely observe the following safety requirements for the connection of an external warning lamp:

To avoid the risk of an electric shock with consequential serious health problems as well as fire hazard of the transformer, if you use an external transformer to power the lamp, please note:

A safety transformer complying with IEC/EN61558-2-6 must be selected.

The internal relay of the device is closed, if the device is in the **READY** operating state

The prescribed laser warning lamp may, of course, also be turned on and off via a normal light switch if the laser is operated.

4.2.5 Foot switch

Now, connect the foot switch to the foot switch connector provided at the rear panel of the device. Make sure to insert the plug as deep as it will go. Then, lock it by turning the cap nut clockwise as far as it will go. The foot switch should be positioned right next to the device.

To be able to activate the foot switch later, the cover must be pressed first so that it opens. Thus, the real switch (black) can be accessed and can be activated by the foot to release the laser emission in the **READY** operating state.


Figure 10: Standard foot switch

4.2.6 Optical fibres

The fibres for endoscopic and surgical procedures (Bare fibres with different fibre core diameters) as well as the transfer fibres for the handpiece and the scanner can be screwed to the device via a SMA connector (see figure below) at the rear side of the device. After that the connector for the fibre identification is to be connected to the device. Via the fibre identification it is par example ensured that the maximum power for fibres with small core diameters is limited.

4.2.7 Focusing handpiece (for QuadroStarPRO^{YELLOW} and QuadroStarPRO^{GREEN})


Figure 11: focusing handpiece

The size of the laser spot can be modified by the use of different optics.

The standard handpiece has a spot size of 1mm. The parts shown in the figure above can be plugged together very easily (spacer to the spot size optics and those to the optics adapter). It is not necessary to screw.

If a different spot size is required, only the spot size optics has to be replaced, the spacer remains the same.

4.2.8 Scanner with integrated skin cooling (optional for QuadroStarPRO YELLOW / GREEN)

The scanner can be easily plugged to the optics adapter of the transfer fibre.

The cable for the skin cooling water can be connected to scanner connector of the rear side of the device (electrical plug and water connectors).


Figure 12: Scanner with integrated skin cooling

4.2.9 Mains connection

Important information

When selecting the power outlet, consider all technical data of your local network and the device. First, verify that the available line power corresponds with the specified power requirements (line power indicated on type plate).

If the device is connected to an inappropriate power supply the device can be damaged.

Please observe the requirements stipulated in your national regulations/standards as amended.

Before the device is connected to the line power supply, the rocker switch of the main switch must be in position O (device without potential). First, plug the mains cable into the power inlet socket at the rear panel of the device before connecting the other end with the mains plug. The operation of the laser requires that the device is connected to a single-phase power outlet protected by a separate fuse rated for at least 1500 VA (circuit breaker type C).

4.3 Switching On

Verify that the following requirements have been met before switching on the device:

- The device is in perfect condition (no damages). No objects are lying at the glass surface of the device.
- The voltage specified on the type label agrees with the line voltage available at the place of installation.
- The power plug of the device has been plugged into an appropriate power outlet.
- The foot switch plug has been properly connected to the foot switch socket at the rear panel of the device and has been screwed down.
- The plug of the door contact is properly connected and screwed to its terminal at the rear panel of the device.
- Fibre is connected properly at the rear side of the device and screwed (SMA and electrical plug). If used, handpiece / scanner are properly connected to the optics adapter of the transfer fibre.
- Fibre / handpiece / scanner is in good order, clean and disinfected – and sterilized if necessary.
- Warning lamps are switched on at every laser room entrance. Doors are closed.

Make sure that all safety precautions have been taken.

Then, proceed as follows:

1. Switch the power switch at the rear panel of the device to I (On).
2. Insert the key into the key switch and turn it clockwise as far as it will go.

While the device is automatically testing essential and safety-relevant parts, the start screen appears on the display. This process will take approximately half a minute.

Then the display for setting of parameters appears as described hereinafter. The device is now ready for operation.

5 Operation of the device

The laser device emits high levels of radiation in the visible and the infrared range. The use of controls or adjustments or performance of procedures other than those specified herein may result in hazardous radiation exposure, in particular in serious damage to the eye.

5.1 General information on operation

The treatment parameters are to be set via the touch screen on the display.


Figure. 13: General design of the display

The modification is very easy and self-explanatory:

- In order to modify a parameter / mode, touch the corresponding button  on the left side of the display. The selected parameter will be framed.

- By touching the buttons  or  on the right side of the display you can increase or decrease the value of the selected parameter.

These settings are normally done in STANDBY. The parameter changes the faster the longer the button is touched – within the technical possible limits, which depend on the parameters.

An adaption (slight change) of parameters is possible in READY also – but stepwise only for safety reasons. That means for each change step the button of the display or of the optional foot switch has to be touched or pressed again.

- Each successful actuation of a button is acknowledged by a short acoustic signal.
- By touching another button (e.g. STANDBY / READY also) the change is saved.

- By touching the **STANDBY/ READY** button  on the display the foot switch is activated. After ca. 2 seconds the yellow laser warning symbol is shown in the low right corner of the display. Up to that point the foot switch must not yet be pressed! Subsequently the emission can be started by pressing the foot switch.

Overview of important buttons of the device

Function	Button	Remarks
Selection		By pressing this button at the left side of the display, the corresponding parameter is selected. You can now modify the selected parameter or switch to another mode.
Increase		The selected parameter can be increased or changed by pressing this button.
Decrease		The selected parameter can be decreased or changed by pressing this button.
STANDBY/ READY		In the READY operating state the laser emission with the selected parameters can be started by activating the foot switch. During the READY operating state the yellow laser warning symbol is shown in the lower right corner of the display. The footswitch is deactivated when the device is in the STANDBY operating state.
TOOLS		General settings of the device can be selected and/ or modified (e.g. pilot laser, optional foot switch, scanner, skin cooling).

5.2 TOOLS menu

By touching the **TOOLS** button in the main menu you can go to the **TOOLS** menu – return by touching the **TOOLS** button again.

The display will show you the following selectable functions of the **TOOLS** menu, which can be selected and modified:

Functions/ Parameters	Display	Remarks
Adjusting the pilot laser	PILOT	The intensity of the pilot laser can be gradually modified.
Adjusting the skin cooling		The intensity of the skin cooling can be gradually modified.
Adjusting the optional foot switch		Setting, which parameter can be adjusted <u>stepwise</u> in READY by the optional foot switch: the INTENSITY (fluence power depending on mode) or the SPEED (repetition rate) or OFF (nothing).
De- and activating of the scanner	SCAN	The scanner can be de-(OFF) or activated (ON). When the scanner is switched ON the parameters for the scanning can be set in the SCAN mode after return to the main menu (SCAN mode) by touching the TOOLS button again – after switching the scanner to OFF the last used mode appears at the display.

Functions/ Parameters	Display	Remarks
energy number of pulses	xxxxxxx yyyyyy	The total energy and number of pulses for all modes are displayed – not resettable.

5.3 Adjusting the modes

By pressing the lowest line on the display the user can select one of several different operation modes. The following modes can be activated: BASIC, EXPERT, BURST, SCAN, CW and VARIX.

The intensive laser emission in the visible and the infrared range of the device can cause thermal changes of tissue structures.

If the focussing handpiece is used in the modes BASIC, EXPERT or BURST, it is absolutely necessary to make sure that **the spot size set / displayed at the display is identically with the spot size used at the handpiece** – otherwise the displayed parameters are wrong!

Tip

If fibers are used (optical fibres without focusing handpiece, for surgical applications mainly) the spot size (core diameter of the fibre) is displayed.

Spot size, Fluence and pulse duration of the laser pulse are related to each other physically. Therefore it is technically in principle not possible, e.g. to realize the largest spot with highest fluence and shortest pulse duration at the same time.

The device 'remembers' the last used settings in each mode and displays the same if one returns into this mode.

5.3.1 BASIC mode

The use of the BASIC mode is very easy. The settings in this mode are suitable for a treatment **using the handpiece and a spot size of 1mm (standard)**. The spot size of 1mm is displayed **only** in BASIC and cannot be changed. Then the fluence and the repetition rate (pulse frequency) has to be set only. Depending on the different fluence value, the device automatically uses the shortest suitable pulse duration when emitting laser radiation, as this is the most common setting for dermatological and aesthetical procedures.

Please use the EXPERT mode, if a spot other than 1mm or pulse durations longer than the shortest possible for the set fluence should be used.

This mode BASIC can be activated **only** if the **transfer fibre** with focusing handpiece is connected.

Function/ parameter	Display	Remarks
Display of the spot size	SPOT	The <u>spot size of 1mm is displayed only</u> but cannot be adjusted.
Adjusting the fluence		The energy density (fluence) can be adjusted and is displayed in J/cm².
Adjusting the pulse frequency		The pulse frequency can be adjusted and is displayed in Hz (pulses per second).
Switching between the modes	BASIC	The different modes will be displayed successively.


5.3.2 EXPERT mode

The EXPERT mode is suitable for all procedures with a different spot size and the same hand-piece and procedures where the pulse duration has to be adjusted manually.

Function/parameter	Display	Remarks
Adjusting the spot size	SPOT	The spot size currently used with the selected hand-piece has to be adjusted correctly.
Adjusting the fluence		The energy density (fluence) can be adjusted and is displayed in J/cm ² .
Adjusting the pulse duration		The pulse duration can be adjusted and is displayed in ms.
Adjusting the pulse frequency		The pulse frequency can be adjusted and is displayed in Hz (pulses per second).
Switching between the modes	EXPERT	The different modes will be displayed successively.

5.3.3 BURST mode

In the BURST mode the pulse are emitted with an alternating sequence of 3ms emission and 3ms pause. The number of pulses depends on the adjusted overall duration. Therefore the settings cannot be adjusted constantly.

Function/parameter	Display	Remarks
Adjusting the spot size	SPOT	The spot size currently used with the selected hand-piece has to be adjusted correctly.
Adjusting the fluence		The total fluence in one burst can be adjusted and is displayed in J/cm ² .
Adjusting the pulse duration		The total pulse duration (emission time) in one burst can be adjusted and is displayed in milliseconds.
Adjusting the pulse frequency		The frequency of the bursts can be adjusted and is displayed in Hz (pulses per second).
Switching between the modes	BURST	The different modes will be displayed successively.

5.3.4 SCAN mode

The scanner is switched ON and OFF in the *TOOLS* mode. The spot size is fixed to 1mm. After switching OFF the device activates to the last used mode.

Function/parameter	Display	Remarks
Adjusting the fluence		The energy density (fluence) can be adjusted and is displayed in J/cm ² .
Adjusting the pulse duration		The pulse duration can be adjusted and is displayed in milliseconds.
Adjusting the intervals between the scans		The scan interval (time needed to set the scanner to the next area) can be adjusted and is displayed in seconds.

Adjusting the distance
between the scan spots

The density / distance between the spots in one scan can be adjusted gradually (highest density 100%, spot by spot).

Adjusting the size and
form of the scan area

FORM

The possible forms and sizes can be selected successively. A small spot size of 1mm can also be selected, in order to avoid changing to the handpiece.

5.2.5 CW mode

In the CW mode the laser emits radiation as long as the foot switch is activated (in READY).

Function/ parameter	Display	Remarks
Adjusting the power		The power of the laser emission can be adjusted and is displayed in W.
Resetting the displayed energy and emission time to 0		The used total energy and pulse duration since the start of the device or the last resetting is displayed and can be reset to 0 (automatically after device switching OFF).
Switching between the modes	CW	The different modes will be displayed successively.

5.2.6 VARIX mode

The VARIX mode is preferably chosen when long pulses are used and display of the energy and emission time used during the treatment is necessary. This is the case for treatment of e.g. varices (endovascular procedures).

Function/ parameter	Display	Remarks
Adjusting the Power		The power of the laser emission can be adjusted and is displayed in W.
Adjusting the pulse duration		The pulse duration can be adjusted and is displayed in milliseconds to seconds.
Adjusting the interval between the pulses		The interval between two pulses can be adjusted and is displayed in seconds.
Resetting the displayed energy and emission time to 0		The used total energy and pulse duration since the start of the device or the last resetting is displayed and can be reset to 0 (automatically after device switching OFF).
Switching between the modes	VARIX	The different modes will be displayed successively.

5.4 Starting the treatment

Before you start a treatment verify again, that:

- the handpiece or the scanner is connected properly to the device.
- all precautionary measures for the laser treatment have been taken.
- if you use the remote door interlock, please check the proper function of this mechanism: Open the door and check that the laser cannot be switched to the READY operating state (an error message will be displayed and has to be confirmed).
- the displayed parameters are adjusted correctly.

The direct and scattered laser radiation can cause serious eye injuries.

During the active laser operation all persons present in the laser room have to wear protective eye goggles or client eye shields appropriate for this laser as specified in chapter *Technical Data*. Other types of laser protective eye-wear do possibly not provide the degree of protection needed!

Before pressing the READY key the attending user must check this.

In the *READY* mode you can unintentionally activate a laser pulse by pressing the foot switch! Therefore, the laser device must be switched to *STANDBY* in treatment pauses and at the end of the treatment.

The emission of the laser beam is indicated acoustically by a warning signal and optically by the flashing emission indicator at the top of the display.

The laser beam may only be released if the handpiece or scanner or fibre is precisely directed onto the skin area to be treated.

In case of dangerous malfunctions the red Laser emergency STOP at the front of the device must be activated immediately. If the emergency problem has been solved, restore the functionality of the device by slightly turning the Laser STOP switch manually anticlockwise to its normal position.

By activating the button READY/STANDBY the device can switch to the READY operating mode.

At the lower right corner the yellow laser warn symbol will be displayed:

During the emission of the laser several safety checks are performed automatically.

The laser emission is started via activating the foot switch.

If the footswitch is pressed, the laser radiation is released from the laser tip of the handpiece. The release of the foot switch instantly stops the emission of laser radiation.

Emission of laser radiation is indicated acoustically by a warning sound and optically by the green blinking laser emission indicator on the display.

During the treatment with the scanner it is possible that the temperature of the cooling probe of the handpiece heats up a little bit, particularly for larger areas. If the upper limit of the cooling temperature is reached, the emission of the laser radiation will be interrupted and the NOTE: *Cooldown handpiece* is displayed. The treatment can be continued as soon as the cooling tip is cold enough again.

CAUTION

If during laser procedures the symbol

or

appear in

the Fluence field, the device has detected a deviation of the output values higher or lower than 20 %.

Higher values can cause undesired side effects, whereas lower values can reduce the effect of the treatment.

If this deviation continues inform the Technical Service.

If you have finished the treatment, switch the laser to the safe *STANDBY* operating state by pressing the *STANDBY/READY* key once more. If you intend to interrupt the treatment for some time, switch the laser to the *STANDBY* state, too. In this operating state, it is impossible to release laser emission unintentionally by the footswitch.

After a longer period of inactivity in the *READY* mode, the device automatically switches to the *STANDBY* mode.

5.4.1 Replacement of the handpiece, scanner or fibre

The replacement of any of the application parts should always occur with the device switched off.

The focusing handpiece and the scanner can be used with the versions of the device which emit yellow or green light only (QuadroStarPRO^{GREEN} und QuadroStarPRO^{YELLOW}). Bare Fibres of different types can be used in principle with all device versions. But at small diameters of the fibre core (0.4mm = 400µm or smaller) the highest possible power is limited.

(There is a special handpiece for cosmetic treatments of nails with the infrared device versions, see attachment.)

5.4.2 Switching off the device

Exit the *READY* state and switch to *STANDBY* by pressing the *STANDBY/READY*-key.

To switch off the laser device, turn the key switch anticlockwise as far as it will go. Only then, you should switch off the power switch at the rear panel of the device. By switching the power switch to '0' the device is completely disconnected from the line voltage.

The direct and scattered laser radiation can cause serious damage to the eye.

Never leave the device unattended while it is switched on. Shut down the device completely by switching the power switch.

After you have turned off the device, make sure you have removed the key from the key switch. Keep the key in a safe place to prevent the unauthorized use of the device.

6 Application

6.1 General instructions

This chapter contains detailed information on the operation and application of the laser device. Of course, the instructions given can never replace the expert knowledge for the application of lasers, the careful study of relevant literature, the personal experience gained under supervision and the critical assessment of the current situation.

But these instructions shall be a useful support for beginners and those not working regularly with the laser.

Therefore, we recommend you to study current literature, especially about therapeutic guidelines and to contact physicians using such type of equipment in order to familiarize theoretically as well as practically with the methods of laser treatment before you start to treat clients yourself.

We would like to assist you in making contact with other users. The competent partner of Asclepion Laser Technologies GmbH will readily inform you about other users in detail.

Requirements for the user

The laser emission of the device can cause thermal changes of tissue structures.

Treatment may not be performed by anyone other than a practising physician (referred to as user for simplification hereafter)

Active medical products may only be applied by users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience.

We recommend to test yourself to summarise und understand important information about the interdependency and relationship between the laser beam and tissue, the connection between the application parameters and the principles of laser safety, as given in the following chapters.

In case of the slightest doubt, please consult an experienced colleague and / or an application specialist of your distributor or of Asclepion Laser Technologies before you start treatment.

Users who want to use this device must have been sufficiently informed about the safe and efficient operation of the laser equipment described in this manual before starting to operate it. Persons employed with or authorized by Asclepion Laser Technologies GmbH give instructions that are based on this manual. This is to be reported in the Device Logbook. The Users must get instructions about the following aspects:

- operation and intended use of the device, practical exercises included,
- operating principle and functionality of the device as well as energies used,
- adjustments of all operating elements,
- applications for the use of the device,
- contraindications and side effects,
- warning labels in all operating modes,

Further training requirements vary in the different countries. The user takes responsibility for the compliance of the training with all requirements of the applicable regional laws and regulations.

5.2 General treatment information

CAUTION

The application of laser radiation on the skin can lead to thermal damages of tissue structures – intentionally according to the application or unintentionally.

The physician is responsible for the precise determination of the application, the treatment zone and the treatment parameters as well as for the correct positioning of the handpiece.

During the treatment the temperature of the cooling probe of the handpiece has to be checked by hand as well as the adjusted parameters at the display.

Tip

In a pre-operation discussion, the client should be informed about all risks that may be caused by the laser treatment and the pre- and post-operative care. Moreover, the client should get all information about the mechanism of laser treatment, alternative procedures, the time required for the treatment, the prospects of successful treatment and about the behaviour after the laser treatment.

You are advised to report this conversation and to make the client sign a special consent form for this briefing. It is also recommended to take photos before and after the treatment.

5.3 Skin treatment with the QuadroStarPRO^{YELLOW} and the QuadroStarPRO^{GREEN}

For more than 20 years, green and yellow lasers have been used as the standard lasers for the treatment of vascular, pigmented und unspecific skin changes, due to their efficacy and their low side effects. Therefore the relevant specialist literature¹²³ was included in the following chapters.

CAUTION

Due to the effect of the laser beam side effects, may be observed in the treated area during the treatment.

The degree of these side effects varies from client to client; normally it is very low. The treatment of darker skin may lead to stronger side effects.

Therefore, a test treatment must be performed in any case in a small area as described below to minimize the side effects and to optimize the effect of the treatment simultaneously.

It is important to bear in mind that the side effects can be stronger for the treatment of larger areas than for the treatment of a smaller test area.

Before the clients get a better impression of the treatment procedure and can better decide for or against a treatment before giving their consent.

The contact surface of the handpiece contacts the client's skin during the treatment. There is the risk of microorganisms being transferred.

The selected handpiece and the scanner must be disinfected before each application and cleaned and disinfected after each application session as described in chapter 8 **Cleaning, disinfection and sterilisation.**

6.3.1 Interaction of the laser beam and the tissue

Light energy can be used to intentionally damage components of the skin or the tissue. The principles of the interaction between light and tissue are the following:

- (1) reflection, scattering and absorption of the light in the tissue,
- (2) the effect of the thermal energy, caused by the absorption of the light (using long pulses) and
- (3) the reaction of the skin to the injury.

The wavelength, pulse duration (duration of the emission) and fluence (energy density) of the laser beam are important factors for the biological reaction.

The goal of the dermatological and aesthetical treatments is the removal of unwanted lesions. Ideally the treatment does not include any unwanted side effects. The **selective photothermolysis** method introduced by Anderson and Parish as a treatment concept in 1983 ensures an effective damage of the target and the simultaneous protection of the skin surrounding it.

Therefore, the following aspects are to be observed for a successful treatment with low side effects:

Absorption of the laser beam – wavelength

The different components of the tissue (e.g. water, haemoglobin, melanin and other pigments) absorb light at different wavelengths more or less effectively. Therefore, the selection of the wavelength influences the selection of the suitable laser and the following probability of unwanted side effects. The absorption of the laser beam must be higher in the target (e.g. the blood in the vessel) than the one in the surrounding skin mainly consisting of water. The higher the absorption in the target the greater the heat generated; the lower the absorption (e.g. in the surrounding skin area), the lower the heating.

Both wavelengths, the green and the yellow one, of the *QuadroStarPRO* (GREEN and YELLOW) obtain a high absorption of melanin and blood (see figure below). Yet, both wavelengths also possess a very low absorption of water, which makes them suitable for vascular treatments and those of benign pigmentations in a selective way. But solid lesions are possible to vaporize also.


Figure 14: Absorption coefficients depending on the wavelength of the laser beam

However, you have also to consider the absorption of melanin in the epidermis of dark or tanned skin that is generally not desired (apart from the treatment of benign pigmented lesions) and can lead to more side effects like e.g. hyper- or hypo-pigmentation (brown or white spots).

To minimize these risks, patients with tanned or darker skin (type IV, V or VI) should be excluded from the treatment. Exceptions can be made for some treatments, like the removal of warts, if the patient was informed and accepted the possible side effects.

Fluence (energy density)

The crucial parameters in every treatment with the laser beam is the intensity of the laser beam. This parameter determines the efficacy of the treatment and the occurring side effects. As pulses are principally used for the treatment of skin changes, the intensity is the energy divided by the spot size given in cm^2 (unit of measurement: J/cm^2). While working with the laser, it should be noted, that the spot size is displayed in mm and therefore has to be converted into cm first.

The relationships are as follows:

Energy density [$\text{Joule} / \text{cm}^2$] (energy density = Fluence)	$= \text{energy} [\text{J}] / \text{spot-area} [\text{cm}^2]$
Power [Watt]	$= \text{energy} [\text{J}] / \text{pulse duration} [\text{sec}]$

Spot size


The spot size affects the depth of impact. For this reason a bigger spot size should be selected when treating a bit larger vessels. Yet, with the high absorption of yellow or green light in blood and the high light scattering in the skin the penetration depth is in principle lower compared to infrared light.

In general, the bigger the spot size, the bigger the skin area that is being heated. Therefore when using a bigger spot size the fluence necessary is usually lower.

Pulse duration

The pulse duration of the laser should be shorter or at least not significantly longer than the thermal relaxation time of the target (e.g. vessel). The thermal relaxation time is the time in which the target emits 50% of the heat to the environment [$\text{TRT} \approx d^2 : 2$; d diameter of the target in mm; TRT is measured in seconds].

Thin structures (vessels) emit the heat much more rapidly than thick structures. To achieve the desired temperature in thin targets it is necessary to deliver the required energy as quick as possible. When the pulse duration is too long the heat loss has to be compensated with a higher fluence, which may lead to stronger side effects.

pulse duration
bigger than
thermal relaxation time

pulse duration
smaller than
thermal relaxation time

Figure 15: Thermal Relaxation Time (TRT) and pulse duration

Target	diameter	thermal relaxation time
Melanosoms	ca. 0,5 – 1 μm	ca. 500 nsec (Ellipsoid)
Extended vessels	10 μm – 50 μm	50 μsec - 1 msec
	100 μm	5 msec
	300 μm	45 msec
	500 μm	
	1mm	(500 msec)

Figure 16: Thermal RelaxationTime (TRT) for different diameters of the target (examples)

For an effective treatment the necessary fluence and spot size are selected first. After that the shortest possible pulse duration is set for light skin and thin structures.

When treating dark skin (skin type III) the melanin in the epidermis should not be heated and therefore a longer pulse duration should be selected – no problem for thicker structures.

For the treatment of thin vessels and dark skin a compromise must be made. The pulse duration should be longer than with light skin and thin targets, but shorter than with dark skin and thick targets. It is also acceptable, as the heat development of thin structures is significantly lower than in thick structures.

Cooling the epidermis

To avoid unwanted side-effects, the effective cooling of the epidermis is necessary (particularly important for treatment of bigger areas, like with the scanner).

Therefore the skin cooling system of the *QuadroStarPRO* is integrated in the scanner. The contact of the epidermis with the cold probe is extremely effective. Moreover, the cooling temperature of the skin cooling system is automatically controlled at any point of time.

6.3.2 Classification of skin types

Skin type		
I (‘Irish’ skin type)	extremely fair skin, blond or red hair, freckles, blue or green eyes	always sunburn no tanning
II (frequent in Central Europe)	fair skin, blond hair, freckles blue or green eyes	mostly sunburn slight tanning
III	light brown or brown hair, brown eyes	sometimes sunburn good tanning
IV (Mediterranean type)	light brown skin, dark or black hair, brown eyes	seldom sunburn fast tanning
V (some Indians, South Americans)	dark skin, black hair, dark eyes	seldom sunburn very fast tanning
VI (some African, Indians)	very dark or black skin, black hair, black eyes	never sunburn, very fast tanning

Figure 17: Classification of skin types according to Fitzpatrick

6.3.3 Application parameters

The **modes** BASIC, EXPERT, BURST and SCAN are provided for dermatological treatments mainly.

In the BASIC mode the standard spot size of 1mm is displayed and cannot be changed (ATTENTION: the focusing handpiece has to have the 1mm optics also!). Additionally the shortest possible pulse duration that is suitable for the selected fluence is automatically adjusted. For the said treatments of light skin these settings are the mostly used ones. Therefore the user must only modify the fluence and the pulse frequency.

As another spot size may be required, and especially for the treatment of dark skin longer pulses are used, the modification of the spot size and other parameters can be made in the EXPERT mode.

The BURST mode enables the radiation in short bursts of pulses, i.e. the selected fluence and pulse duration is dispersed in small bursts of pulses. This way, the treatment, especially of darker skin, is gentler.

The very important parameter for the treatment of all skin lesions is the **FLUENCE**. This parameter can be modified on the display in all modes suitable for skin treatment. The higher the fluence, the more effective is the treatment, although stronger side effects are also possible. For this reasons test treatments have to be made and evaluated before the treatment itself!

The **spot size** depends in the used handpiece and is also displayed. In the BASIC and the SCAN mode it cannot be modified and is automatically set to 1mm. This spot size is suitable for most of the treatments. The spot size of 0.5mm is even gentler for the treatment of very small and superficial vessels (shorter pulse duration at the same fluence is possible). With this small spot size the user has to pay special attention to the target, as it may be more difficult to take aim.

When treating deeper lesions, a bigger spot size may be suitable. It should be noted, that treatments using a bigger spot size require a lower fluence and should use a lower fluence also, as a larger area is heated. Moreover the performance of test treatments is strongly recommended.

Special attention has to be paid to the correct setting of the spot size diameter (see marking on the handpiece). An incorrect spot size leads to false adjustments of the fluence and the pulse duration and therefore can lead to overtreatments.

The **pulse duration** can be modified as explained above in the modes EXPERT and BURST. When treating thin vessels and benign pigments it should be selected as short as possible (with a sufficient fluence and spot size). With darker skin a longer duration is possible, in order to treat this skin types more carefully. Short breaks between the bursts (BURST) also lead to a gentler treatment of the melanin in the epidermis.

The selected **pulse repetition rate** in Hz (in BASIC and EXPERT mode) and the pulse duration in ms (in BURST and SCAN) determine the duration of the treatment. Both have no direct impact on the efficacy of the treatment. In very sensitive areas a reduction of the treatment velocity may reduce any pain.

The **Scanner** is suitable for treatment of bigger areas. It treats a selected pre-selected **FORM** of an area automatically. The distance between the spots in the treatment area can (**density**) be modified and adjusted for different treatment techniques: spot-by-spot or polka-dot technique (gaps between spots, see below Application techniques). Beside FORM and density the **fluence** and **pulse duration** as well as the **interval** (time to move the scanner to the next place) can be set.

In general: The higher the fluence, the larger the area, the more dense the spots the more absorption by the lesion itself – the more side effects have to be expected.

A sufficient **skin cooling** may reduce side effects when treating with thermal operating techniques. The treatment areas are often cooled with cooling pads. The bigger the treatment area the higher the heat development and the higher the risk for unwanted side effects. For this reason the scanner, which treats bigger skin areas, is equipped with an integrated skin cooling.

6.3.4 Application techniques

The selection of a suitable application technique is as equally important for the clinical and aesthetic result of the laser treatment as the correct selection of application parameters.

Single spot technique

- Technique for coagulation of single vascular, pigmented or solid lesions (e.g. syringomas)
- Use a suitable spot size for every individual localisation

Single vessel technique

- Use this technique for the treatment of visible single vessels. Follow the course of the vessel continuously without any overlaps.
- Continue the treatment until the required results are met.
- Do not treat the interstices between the vessels.
- If necessary treat the main vessels first, before continuing with the peripheral smaller vessels (e.g. with spider nevi)

Spot-by-spot technique

- Use this technique for smaller vascular or pigmented lesions in skin areas with a high regeneration capacity and a low scar risk.
- No overlaps
- Technique for vaporisation

Polka-dot technique

- Gentlest treatment technique, as the healing starts from the healthy interstices between the spots – the bigger the distance between the spots, the gentler the treatment.
- Use this technique for treatments of strong absorbent changes which tend to stronger side effects (e.g. big vascular lesions) or in sensitive regions with higher risks.
- Close the gaps in the second session. After the first 2 sessions the absorption is usually significantly lower and the reaction of the skin or tissue is known, so that in the following sessions (if necessary) a technique without gaps can be used normally (e.g. spot-by-spot technique).

Scanner

- The scanner is mostly used for the treatment of extensive lesions, but can also be set to a smaller spot size and used like the handpiece.
- Besides the adjustment of the different forms and sizes, the distance between the spots can also be modified, like for the spot-by-spot or the polka-dot technique, depending on the intended treatment.

6.3.5 Selection of patients

WARNING

Unqualified laser treatments may lead to serious health damages.

The professional qualification including practical experiences with the planned application is a crucial requirement for the attending physician and the procedure.

The physician is responsible for the precise determination of the application for the given indications.

The correct selection of patients is a crucial requirement for an appropriate efficacy. The patients have to be informed about the method, risks, alternative procedures, possible unwanted side effects and the prospects of success. Possible preoperational treatments of the lesions and their side effects should also be explained. Any side effects of a pre-treatment have to be healed before any test treatments can be made.

Patients with an unrealistic expectancy should not be treated.

Patients with the skin type I to III, who are **not tanned** in the area to be treated, can be treated with the following restrictions.

- When treating patients with skin type III the parameters for the test treatment should be adjusted very carefully (long pulse duration and low fluence).
- The indications with the best results with a selective laser treatment are red and extended vessels up to a diameter of 1mm, epidermal pigmentations and other skin changes as specified in chapter 1.3 Intended use.
- It should be noted, that the depth of impact of green and yellow light is restricted to superficial treatments, due to the high absorption in blood and melanin and the relatively high scattering in the tissue. Therefore some lesions like thick port wine stains, haemangioma, bigger vessels and thick or deep pigmentations cannot be treated successfully or completely.
When treating extended vessels (only red and superficial up to 1mm!) on the legs the prospects of success are lower and the patients have to be informed thereon.
- For the treatment of the most of vascular or pigmented lesions, several treatments are necessary (see chapter *Follow up treatment*). Test treatments are always recommended (see chapter *Test treatment*)!

The (unintentional) laser treatment of precancerous or even malignant lesions poses a significant risk for the health of the patients.

It is an absolute prerequisite that in case of the slightest doubt about the benignancy of pigmentation, the diagnosis has to be confirmed or treatment has to be cancelled.

Tip

Heavy smoking and/ or adiposity may influence the result of a vessel treatment.

- Patients with big or suppurating haemangioma in health threatening areas (e.g. mouth, throat, eye) or near the eye and haemangioma that penetrate deep into the skin should be immediately send to a specialised hospital with experiences with these conditions.
- When treating children, the attending physician has to decide whether the treatment is necessary or not, with regards to ethical issues. Especially for infantile haemangioma a profound and an assured diagnosis is indispensable in order to determine a suitable treatment. The treatment parameters and techniques of a child's skin have to be selected very cautious.

- In case children should be treated it is of course to decide by the physician according to ethnics principles, if a treatment is necessary or reasonable.
Especially in hemangiomas of children a well-founded diagnosis is absolutely necessary to determine the most appropriate therapy.
The procedure has to be chosen especially carefully due to the children's skin.

6.3.6 Contraindications

See chapter 1.3 *Intended use*.

6.3.7 Possible side effects

All side effects are less significant for light skin and smaller lesions (e.g. with few and small vessels) [1], [2], [3] and become more intensive for thicker and bigger lesions with darker skin – provided that the skin is not tanned – and with a too high fluence value in the individual treatment (test treatment!). **Skin cooling after treatment reduces side effects!!**

Therefore the start parameters for the test treatment with patients with skin type III should be chosen lower.

During the treatment the clients have a pricking feeling or even burning sensation connected with the heat when treating thick lesions or extensive skin areas.

After the treatment erythema and oedema appear in the treatment area, e.g. on the course of vessels. The skin around the eyes and on the neck is very sensible [2]. These effects are a normal reaction of the skin and the tissue and usually last a few hours or days. After some treatments of haemangioma these effects may persist up to a few weeks.

Increasing (instead of healing) swelling, reddening, crusts, pain and a high temperature are symptoms of an infection. Topical antiseptics or oral antibiotics should be subscript [2].

Blisters, crusts and transient pigmentary changes (last up to 6 months) occur occasionally but can be reduced to a minimum by test treatments, sufficient skin cooling, protection from light after the treatment and by not treating tanned skin [1], [2].

Hyperpigmentation often occurs with darker skin but can be reduced or avoided with a pre-treatment with bleaching cream (e.g. treatment with hydrochinon for two weeks). Hypopigmentations are usually caused too high parameters and may persist, especially near the neck, on the legs and the chest (test treatment!) [2].

After the treatment of haemangioma crusting is usually and persists up to 8-10 days.

In treating leg veins hyperpigmentations occur more often and can persist for months [1]. With the yellow or green light one should treat small red veins only; in this case the risk is lower compared to larger veins.

The reactivation of herpes simplex (especially with patients with herpes recurrences more often than 6 times a year) when treating areas in the face or on the legs (in case of Herpes simplex genitalis) should be prevented by subscribing oral virostatic agents one day before the treatment [2].

Blisters and crusts should be treated prophylactically with antiseptics [3] and usually heal within a week and leaving now scars. Pigmentary changes may persist for a few months.

Changes of the skin texture and sunken scars (usually caused by too high parameters) are very seldom for the most of the indications with appropriate procedures and patients' compliance [1], [2], [3]. However, these side effects may occur after the vaporisation and coagulation of lesions, like Xanthelasma, haemangioma and especially warts more often [1]. The patients have to be informed hereon. Usually the patients accept these side effects, as there is still a significant improvement.

The development of granuloma pyogenicum has been reported in some cases [3]. There is a low risk of a superinfection [1].

Nonresponders are possible for some indications, e.g. lesions with a high blood flow. Patients have to be informed accordingly.

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[2] Vascular lasers and IPLs: Guidelines for care from the European Society for Laser Dermatology (ESLD)

M.Adamic et al., Journal of Cosmetic and Laser Therapy, 2007;9:113-124

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6.3.8 Preparing laser treatment

As a compulsory pre-condition, the area to be treated must be absolutely shaved and cleaned. Remove also deodorants or other cosmetic products.

To document the treatment success photos should be taken.

Anaesthesia is usually not necessary (only with bigger lesions) if the skin is cooled sufficiently. When treating smaller vessels a local anaesthesia should be avoided, as the absorbent blood is reduced by the vasoconstriction and therefore the efficacy of the laser treatment might be reduced.

Suitable safety goggles, eye shields for the patient and cold cool packs (from the fridge, must not be frozen) should be provided. A surgical mask and a suction device are necessary when vaporising tissue, e.g. in treating and warts.

The user should mark the area to be treated with a white kajal, as it is difficult to localise the vessels when erythema occur during the treatment.

Pre-cooling the treatment area (e.g. with cool packs or the cooling area of the scanner) is very effective and therefore recommended for the reduction of side effects and discomfort, except when treating epidermal pigments.

6.3.9 Test treatment

The test treatment should always be performed with the same application parameters and techniques including cooling before and after the procedure as the treatment itself.

If the patient decides in favour of the laser treatment after the medical briefing, a test treatment should always be performed so as to determine the optimal parameters for the individual situation and to evaluate the efficacy and the possible side effects. In order to eliminate of the development of pigmentary changes the user should wait an appropriate time before evaluating the side effects.

When using the scanner for the test treatment the spot size can be set to a small area.

Before beginning the test treatment a suitable spot size (usually 1mm) is adjusted. For the treatment of very small vessels a smaller spot size might be chosen and for treatments of lesions that deeper in the skin, a bigger spot size is suitable.

When treating bigger and thicker haemangioma (max. thickness ca. 1-2mm) it is gentler and therefore strongly recommended to treat through a flat ice cube with no bubbles or a very cold glass plate. In pressing the handpiece on the ice cube or glass plate the distance to the skin can be reduced and the skin is cooled simultaneously. The scanner has an integrated skin cooling, which also should be used with a light pressure. The same should be applied when treating papular spider veins as they are more difficult to treat and longer pulse duration should be used.

When treating patients with light skin and using the standard spot size of 1mm, the BASIC mode should be chosen (this automatically adjusts the shortest pulse duration for the given fluence). For patients with the skin type III the EXPERT mode should be chosen, as this enables the selection of longer pulses and therefore a gentler treatment of darker skin. When using another spot size as 1mm the EXPERT should also be used and the according spot size should be adjusted correctly.

In the SCAN mode the same parameters as in the EXPERT mode can be adjusted and additionally the desired spot size, form and density. The most important additional setting is the distance between the different spots in the treatment area. Similar to the polka-dot technique the distance between the spots can be chosen to be bigger. This provides a gentler treatment, which is very important when treating with the scanner.

CAUTION

The laser radiation can cause unwanted side effects, especially when used with incorrectly adjusted treatment parameters. Therefore special attention has to be paid when adjusting the correct spot size.

If the spot size is not set correctly, the adjusted fluence and pulse duration are incorrect!!!

If the spot size is set too big, the value of the fluence will be too high (risk of stronger side effects). If the spot size is set too small, the value of the fluence will be too low (no result).

When the test treatment is prepared properly the test procedure should be started in an insensitive area with low fluence.

While performing the test treatment the user should pay special attention to ensure that:

- the warning lamp of the door is on,
- all people in the laser room wear appropriate safety goggles and
- the eyes of the patient are protected appropriately.

Laser radiation can cause serious eye injuries.

Everyone present in the laser room during a treatment session must wear laser safety goggles including eye protection for the client according to the instructions given in chapter *Technical Data*

The set fluence as well as the temperature of the skin cooling probe must be controlled before and during the treatment!!

Then the user should proceed as follows:

- switch the device to READY,
- check again the displayed parameters,
- place the handpiece or the scanner on the skin to be treated,
- activate the foot switch for a few spots (low pulse frequency)
- evaluate the result.
- If no result is visible, the fluence can gradually be adjusted until the desired effect is apparent or
- If the effect is too strong, reduce the fluence.

Clinical endpoint:

- Ectatic vessels: Blanching (disappearing of the vessels) of greyish blue discolouration - a persisting white discolouration (scar) caused by a too high fluence value must be avoided!
- the greyish blue discolouration of the vessel (coagulated blood is darker) that can occur with bigger vessels may persist for several months – the lymphatic system needs time!
Note: Different diameters of a vessel cause a different blood flow and therefore a different absorption. For this reason the test treatment has to be repeated when the diameter changes
- Epidermal pigmentations: greyish brown discolouration (debris), which can be removed with a saline solution moistened pad simply before the evaluation of the effect

- Xanthelasma: greyish white discolouration, can be removed with a saline solution moistened pad
- Warts and similar lesions: carbonisation, if a vaporisation is wanted (virus protection!), greyish white discolouration in case of coagulation (peeled wart)

The treatment areas have to be cooled after the test treatment.

The device has to be switched to STANDBY after the treatment.

A test treatment can help to evaluate the risk of pigmentary changes and of cicatrisation as well as the healing process. Additionally the patient can also see how the lesions are treated.

When treating pigmentations and patients with darker skin the evaluation has to follow after at least 4 weeks, in order to assess the possibility of pigmentary changes.

6.3.10 Treatment procedure

Once the treatment area has been prepared and the parameters are adjusted according to the results in the test treatment, the lesions can be treated. Depending on the size of the lesions single spots are applied or the area is treated with the spot-by-spot technique (without overlaps). When a very gentle treatment is required (e.g. with port-wine stains) the polka-dot technique is recommended (can also be used with the scanner).

When treating single vessels the single vessel technique is used (no overlaps and no treatment of the interstices between the vessels)

When treating spider nevi the main vessel should be treated before continuing with the smaller vessels. A larger spot size might be necessary for the treatment of bigger vessels.

When the treatment of one area is finished a cooling pad is applied before continuing with the next area with the same technique. The bigger the treatment area, the more heat is applied and the stronger the reactions of the skin. Therefore treatment breaks should be taken to supply sufficient cooling of the skin or if necessary the treatment should be spread over several sessions.

The device has to be switched to STANDBY after the treatment.

The whole treated area has to be cooled for at least 10 minutes after the treatment, depending on the size of the treated lesions. This proceeding and the continuing cooling ensure a significant reduction of the sensation of pain and the side effects.

6.3.11 Post treatment

After the skin has been cooled sufficiently, a skin care product or if necessary an antiphlogistic can be applied. Continuously cooling can considerably reduce the side effects, even when treated with slightly too high parameters. The bigger the treated area, the longer the skin should be cooled. For example, you can recommend your patient to cool the area with a cool pack in a soft tissue for 15 min every hour for four hours [2].

If an overtreatment occurred topical antibiotics should be applied and the healing process should be monitored.

CAUTION

If the treated area is exposed to sun and solar radiation after the treatment, stronger side effects can develop, in particular pigmentary changes.

Therefore, such a radiation should be avoided for a period of 6-10 weeks or at least a sun screen with a protection factor not below 50 should be used.

The lased area should not be traumatized either for 2 to 3 days.

When treating near the eye (e.g. on the cheeks) an oedema may occur under the eyes after 1 to 2 days. Recommend your patients to sleep with the head elevated on a cushion.

After treatment of ectatic vessels, patients should refrain from activities that cause the body or the vessels additional physical strain, like jogging, going to the gym, sauna sessions, hot baths and similar activities.

Patients must not manipulate any developing crusts.

After the treatment of haemangioma, the treated area usually dehydrates and develops crusts which persist for 8 to 10 days. Oedema and small erythema may persist up to a few weeks. Thick haemangioma are difficult to be treated as the green or yellow light may not penetrate the skin deep enough.

6.3.12 Follow up treatment

When the polka-dot technique was used in the treatment a second session is necessary to treat the interstice spots. After the second session the polka-dot technique usually is not necessary any more. As the absorption is less strong the side effects after follow up treatments are also less severe.

Usually, the final treatment results are visibly after 4 to 10 weeks and after the treatment of ectatic vessels of the legs it may take 3 to 6 months. You can easily check with finger pressure whether the blood is flowing or is coagulated and therefore is still being reabsorbed. For this reason follow up treatment should be planned after a suitable period of time.

Small vessels in the face can usually be removed successfully with 1 or maximum 3 treatments. Yet some procedure may require several treatment sessions every 4 to 6 weeks in order to heal properly and achieve a satisfying treatment result. For some indications, e.g. port wine stains, a successful removal or bleaching cannot be guaranteed. The patients have to be informed hereon. Extensive port wine stains should be treated with a dye laser.

Ectatic vessels near the lower ala of the nose and in the sebaceous skin near the root of the nose or the chin are often more therapy resistant. Those indications might relapse within only a few weeks or need several treatments.

6.4 Treatment with the infrared version of the *QuadroStarPRO*

For over than 20 years now the Nd:YAG laser (1064nm) is an established standard lasers in the surgical laser therapy besides the CO₂ laser. In the 90s the diode laser became more popular and continuously replaced the Nd:YAG laser, at first with the wavelength on 810nm and later with 940nm and 980nm.

These different wavelengths in the near infrared range (NIR) are all very suitable for surgical procedures for the coagulation and vaporisation of various types of soft tissue. Their mode of action is similar and they can also be used for endoscopic procedures contrary to the CO₂ laser.

Unqualified laser treatments may cause severe health damages.

The professional qualification including practical experiences with the planned application, knowledge about the indications and techniques as well as the suitability of the surgical room and instruments are a crucial requirement for the attending physician and the surgical or endoscopic procedure. The physician is responsible for the sterility of the procedure in which the laser is usually only a small part. Moreover he or she must be able to finish the procedure conventionally in the case of an emergency.

The physician is responsible for the precise determination of the application and the method for the given indications.

The NIR versions of the *QuadroStarPRO* are usually used with an optical fibre (and without a focusing handpiece or scanner). Sometimes a handpiece without optics, but with cannula as a holder.

During the procedure the fibre usually has contact with the tissue of the patient. There is a risk of transmitting microorganisms.

The fibre and cannula has to be sterile before each treatment and in the case of skin treatments at least disinfected.

If the end of the fibre has contact with any kind of metal, the laser beam causes a heating. There is a fire hazard! Please pay attention to the warnings given in chapter *Safety* and especially the following:

- If you use endoscopes, trocars, cannulas or similar devices, always pay attention, that they are not inflammable and that the end of the fibre always protrudes the instrument at least 7 to 8 mm during laser emission!
- For procedures where the fibre is not in constant contact with the tissue, the end of the fibre has to be kept clean!
- For all endoscopic, endovascular, interstitial and similar procedures only fibres with a thermally stable cladding (e.g. quartz-quartz-fibres) may be used!

In case of an emergency immediately activate the Laser-emergency STOP switch in order to immediately stop laser emission.

6.4.1 Interaction of the NIR laser beam and the tissue

The general connections and interactions are the same like the ones for skin treatments explained in chapter *Interaction of the laser beam and the tissue*. Yet, the mode of action due to the different wavelength and therefore the possible procedures and techniques are different.

Absorption of the laser beam – wavelength

The effect of NIR laser radiation to the skin or tissue is different to that of yellow or green wavelengths. The absorption of the light in the blood and the melanin is significantly higher and in water lower. Additionally the scattering of the light in the tissue decreases with an increasing wavelength (see Absorption curve in the chapter above).

These characteristics allow the NIR laser radiation to penetrate deep into the tissue and, unlike the yellow and green wavelengths, enable a deep coagulation and vaporisation of different types of tissue. Depending on the energy- and power density the tissue can be coagulated or vaporised.

Using the fibre contact procedures it is possible to cut tissue without causing a bleeding. The diode laser with the NIR radiation is also very suitable for the endovascular coagulation of varicose veins.

Due to the decreased selectivity compared to the yellow and green versions of the *Quadro-StarPRO* (blood and melanin absorption significantly lower) the NIR laser is not suitable for the treatment of small vessels and epidermal pigmentations. On the other side, the risk of pigmentary changes is significantly lower after skin treatment compared to green / yellow laser light.

The effects of all NIR lasers with the wavelengths from 800 to 1100nm are in general very similar. Due to the slightly higher absorption in haemoglobin a laser the wavelength of 940nm is usually preferred, e.g. for recurrent nose bleeding or endovascular procedures.

The absorption of water is the highest at 980 nm wavelength. The absorption by haemoglobin is also slightly higher than with the Nd:YAG laser. Therefore this wavelength is very suitable for coagulation, vaporisation and cutting of soft tissue in surgery. Due to the heating skin firming effect can also be noticed, e.g. with laser assisted liposuction.

In the dental field most recently a laser with the wavelength of 980nm is used for the treatment of soft tissue as well as for tooth bleaching.

Power – power density

For the treatment of skin and in surgery the crucial parameter is the intensity of the laser beam. This parameter determines the efficacy and the following side effects. In contrast to skin treatments, the NIR laser beam is usually radiated continuously or in relatively long pulses and a focussing handpiece with a set spot size is usually not used, the suitable parameter to be adjusted is the power of the laser in these cases.

Analogous to the fluence during skin treatments the intensity of the laser beam is the power density. That means that the size of the irradiated area is important because it is squarely included.

The mathematical relationships are as follows:

$\text{Power density (Watt/ cm}^2\text{)} = \text{power [W] / irradiated area [cm}^2\text{]}$

The irradiated area depends on the distance between the tissue and the fibre and can be recognised by the size of the pilot laser spot.

The shorter the distance of the end of the fibre to the tissue, the smaller is the irradiated area. When the fibre is in contact with the tissue the area is the smallest size and is as big as the fibre's diameter. The thinner the fibre, the higher is the power density at the same power

Therefore the intensity of the laser beam can be modified by changing the distance between the fibre and the tissue, e.g. coagulating with a bigger distance, vaporising with a smaller distance or even cutting when in contact.

For very delicate and small the 400µm fibre is better suited than the 600µm fibre; less power is used.

When working with a non contact technique the end of the fibre must stay clean all the time. Any small pollution will absorb the laser beam and cause a heating or it even glows at the end of the fibre. When working with a contact technique the cutting only works properly when the end of the fibre is blackened. The end of the fibre is cooled by the constant contact to the tissue.

Duration of the radiation

When using the laser beam for surgical purposes it has to be noted absolutely, that the radiation dehydrates the tissue up to a carbonisation. The blackening of the tissue abruptly increases the absorption and therefore the effect, e.g. the coagulation can be very slow, but when the tissue is blackening a perforation can be initiated very quickly (on purpose or unintentionally).

This general rule can be applied: The lower the power density and the longer the radiation without carbonisation, the bigger is the coagulation zone (the deeper the necrosis). When vaporising the coagulation zone is limited and depends on the power density and the duration of the radiation.

If the necrosis has to be small, a short distance or contact method should be chosen and the radiation should be emitted in pulses or the fibre is moved quickly over the tissue to be treated.

Types of fibres

The core of the optical fibre is made of quartz of different diameters. Usually a diameter of 600µm is used. When cutting with a contact method (higher power density and smaller necrosis) or for example with an endoscope with the Albarran lever (more flexible) even smaller diameters are possible. The maximum possible power with diameters smaller than 400µm is limited, but less power is necessary for a high power density.

The cladding of the fibres can be made of different materials, which may have different apertures and characteristics. Because of safety reason we deliver Quartz-Quartz fibres only (no danger that the plastics of the cladding „burns back“).

The numerical aperture is the angle of the fibre's opening. The smaller the aperture the lesser the size of the irradiated area changes when using a different distance.

Example:


Intensity	= Power density
	- Important parameter for effect on the tissue
Power density [Watt / cm ²] = power per irradiated area (Spot)	
For example:	Power 10W
	Fibre 600µm - Contact (0,003cm ²): 3600W/cm ²
	Fibre 400µm - Contact (0,001cm ²): 10000W/cm ²
Vice versa	In order to achieve 10000W/cm ² you need (with contact):
	At 400µm- Fibre: 10W
	At 600µm- Fibre: 30W
	
	PCS Fibre: bigger aperture (More sensitive to changing the distance)
	600µm - Distance 1mm: Spot 1,36mm - 690W/cm ²
	600µm - Distance 5mm: Spot 3,60mm - 98W/cm ²
	Quartz - Quartz - Fibre: smaller aperture
	600µm - Distance 1mm: Spot 0,96mm - 1383W/cm ²
	600µm - Distance 5mm: Spot 2,40mm - 222W/cm ²

Figure 18: power density – types of fibres

Tissue and application environments

The diode lasers can be used in different medical specialties and with different types of tissue. As the characteristics of the tissues can be very different, the parameters for the same treatment may also vary and need to be adjusted specifically.

Depending on the application the fibres may be used in different environments, like air or liquid environment (e.g. irrigation treatments in the urology). When used in a liquid environment the cooling effect of the liquid on the fibre and the tissue has to be considered when adjusting the parameters.

6.4.2 Application parameters

The modes CW and VARIX are provided mainly for surgical procedures.

The CW mode is the easiest to use and is usually used for all standard procedures. Only the power value has to be adjusted and the laser emits radiation as long as the footswitch is activated (in the READY state).

In the VARIX mode the power, the pulse duration and the interval between the pulses can be adjusted. This mode is suitable for procedures with short interruptions, like the endovascular treatment of varicose veins. During the pulse interruptions (intervals) the fibre is pulled back.

This mode may also be used to adjust the pulses manually. Additionally this mode displays the total applied energy and emission time, which can also be reset to 0 (if a new patient has to be treated).

For surgical treatments the Power is a crucial parameter. The other parameters that influence the effect of the procedure are explained in the chapter before.

The pulse duration (and interval) used in surgical procedures influences the necrosis zone necessary (e.g. cutting with contact method and with short pulses causes a small necrosis) and the treatment procedure as it is explained with endovascular treatments above.

NIR laser radiation can cause severe health damages, if certain structures are damaged unintendedly.

The parameters should be chosen as low as possible for the desired effects from experience. Due to the different influences, the desired tissue alteration is the key indicator for the right choice of parameters.

In particular it has to be remembered how deep the coagulation is. It has to be absolutely excluded that an unintended perforation is caused during healing process due to resorbed tissue !!

6.4.3 Application techniques

Soft tissue can be **cut** by the means of the **contact technique** and by placing and pulling the fibre over the tissue. No pressure is needed. The size of the necrosis zone depends on the fibre's diameter (400µm fibres have a smaller necrosis zone than 600µm fibres), the adjusted power and the speed with which the fibre is pulled. Usually the necrosis zone is only less than a millimetres small. The laser should be used in the CW mode and for working with smaller structures and fibres lower values are to be adjusted or pulses are to be used (VARIX).

When a thicker necrosis zone is needed for tissue with a high blood circulation or tissue should be removed a **vaporisation** with a non-contact method can be used. A smoke evacuation system is necessary for this technique. When vaporising, a higher power density is needed than the one applied for coagulation. That means a smaller irradiated area (smaller distance), a longer radiation time and / or a higher power value should be used.

It should be noted, that the laser effect develops relatively slowly as long as the tissue has a light colour. During the coagulation, the tissue dehydrates and carbonizes. In the moment of the carbonisation the absorption changes abruptly because of the blackening and the following vaporisation is much faster. Usually the necrosis zone around vaporised tissue is bigger than around cut tissue and smaller than coagulated tissue.

Usually a non contact technique is used for a deeper **coagulation**. If the fibre comes in contact with the tissue the end of the fibre is polluted and needs to be cleaned before the non contact treatment can be continued. The end of the fibre must be kept clean at all times! Otherwise the laser light will be absorbed by the particles on the end of the fibre. As no tissue contact and therefore no cooling is provided this may lead to a heating or even annealing of the fibre. Additionally only a part of the laser beam will reach the tissue.

The distance usually is 5 to 10mm and the standard fibre is the 600µm one. The bigger or deeper the coagulation should be, the bigger the distance and the longer the radiation time. A vaporisation should be avoided as the coagulation near the vaporised tissue is limited.

Another possible technique is the **interstitial coagulation**. Usually only low power values are necessary, as a carbonisation should be avoided. Occasionally, after the pricking with the fibre (e.g. by the means of cannulas or similar instruments), small spherical coagulation necrosis are placed around the end of the fibre (e.g. when treating liver metastases).

The fibre also can be pulled back after the pricking while emitting laser radiation. If a bigger area is treated this can also be repeated several times, e.g. like a fan-shaped method for the laser assisted liposuction.

In some procedures (e.g. thick haemangioma) the pilot laser is visible through the skin and the location of the fibre end can easily be controlled, otherwise an ultrasound device should be used to ensure the localisation of the fibre end (e.g. with endovascular treatments).

In order to fully understand and visualised the general mode of action the procedure may be tested on poultry meat, beef or pork including fat and the liver (bright and dark muscle tissue, fat and the liver) using different application parameters and different procedures. If possible material more similar to human tissue should be used.

In this way the effect of the same parameters on different tissue and of different parameters on the same tissue can be tested and demonstrated.

Naturally these tests cannot substitute professional knowledge. A visit with a specialist who is an expert for this technique is strongly recommended in order to provide the surgical laser procedures with the necessary safety.

Examples for different applications

(Summarised and adapted* for the *QuadroStarPRO* from:
Berlien, Müller: *Angewandte Lasermedizin*, ISBN 3-609-70510-8.)

* due to max. power 25 W the size of the lesion to be treated is limited

Effect	Application	Parameter	Remarks
Coagulation in air small coagulation zone (1 – 2mm) e.g. for: endoscopic procedure for M. Osler or vascular malformation	Bare fibre 200 – 600µm Non contact In air	Power: 25 W* Spot: 2mm Pulse duration: 0,2 – 0,5s Interval: min. 0,2 sec	- Clean and freshly broken fibre - Avoid carbonisation (irrigation, if necessary) - Avoid dehydration - Avoid popcorn effect - repeated exposition if necessary - Stop when tissue begins to blanch
Coagulation in water small coagulation zone (1 – 2mm) e.g. for: bladder polyps, intraurethral Condylomata acuminata	Bare fibre 200 – 600µm Non contact In water	Power: 25 W* Spot: 2mm Pulse duration: 0,2 – 0,5s Interval: min. 0,2 sec	- Clean and freshly broken fibre - repeated exposition if necessary - Stop when tissue begins to blanch
Vaporisation in air with small coagulation zone (0,5 – 2mm) e.g. for: precise endoscopic resection of scars or tumours of the respitorial / digestive tracts, during laparoscopic pro- cedures, especially for children	Bare fibre 200 – 600µm Contact In air	Power: 15 - 25 W* Pulse duration: 0,1 – 0,3s Interval: min. 0,1 sec	- Preblackening of the fibre - Coagulation zone depends on the relation between power and pulse duration - No drilling, place fibre without pressure - Fibre tip has to be visible at all times
Vaporisation in water with small coagulation zone (0,5 – 2mm)	Bare fibre 200 – 600µm Contact In water	Power: 25 W* Pulse duration: 0,1 – 0,3s Interval: min. 0,1 sec	- Preblackening of the fibre - Coagulation zone depends in the relation between power and pulse duration - No drilling, place fibre without pressure - Fibre tip has to be visible at all times

Effect	Application	Parameter	Remarks
Laser induced coagulation Interstitial Coagulation (LIC) e.g. for: Subcutaneous vascular malformations, Haemangioma, Perforating varices, ITT of malignant tumours, Metastases	Bare fibre 600µm contact	Power: 5 W Pulse duration: cw (continuously)	- Clean and freshly broken fibre - Use 16G cannula made of Teflon or metal - Fibre protrusion 8-10mm in order to avoid heating or melting! - Position of fibre visible through pilot laser, ultrasound or MRI - Only retrograde procedure! - Crepitation shows best tissue reaction, palpable through the skin or vibration of fibre - Process monitoring with coloured duplex ultrasound or MRI temp. measurements - Coagulation zone depends on pulse duration
Interstitial/intraluminal coagulation (LIC with irrigation e.g. for: intravascular treatments pathologically changed vessels, e.g. vascular malformation)	Bare fibre 600µm contact	Power: 8 - 10 W Pulse duration: cw	- Irrigation with 100mm H₂O and 2ml/min - Same procedure as described above
Interstitial coagulation (impression technique)	Bare fibre 400 - 600µm Contact In air -	Power: 5 - 8 W Pulse duration: cw	- Freshly broken fibre - Coagulation zone depends on pulse duration - Place end of the fibre with light pressure - Spot sizes coagulation of the tissue surface, volume coagulation under the tissue surface
Fistula shrinking: e.g. for fistula in the respiratory or digestive tracts, including post operational, congenital and acquired (e.g. anal fistula)	Bare fibre 600µm	Power: 4 - 5 W Pulse duration: cw Retraction speed: 0,5 - 1mm/sec	- Freshly broken fibre - Pilot laser shows fibre position - Only retrograde procedures - Additionally fibrin glue if diameter more than 2mm - Maximum diameter 5mm - Avoid vaporisation

Effect	Application	Parameter	Remarks
Coagulation, slight bleeding bigger coagulation zone (2 – 3mm)	Bare fibre 400 – 600µm Non contact In air	Power: 25 W* Spot: 2 - 5mm Pulse duration: 1 – 2sec. per pulse Retraction speed: 0,5 - 1mm/sec	- Freshly broken fibre - Avoid carbonisation - Irrigate with sodium chloride in case of bleeding - repeated exposition if necessary
Coagulation, stronger bleeding bigger coagulation zone (3 – 5mm) e.g. for: ulcer bleeding, tumour, vascular malformation	Bare fibre 600µm Non-contact In air	Power: 25 W* Spot: 2 - 5mm Pulse duration: cw	- Clean, freshly broken fibre - Use additional gas cooling - Avoid carbonisation - Continuing irrigation with sodium chloride during procedure - continuing exposition for several seconds - avoid popcorn effect, if necessary stop treatment

6.4.4 Selection of the patients, contraindications, side effects

The infrared radiation can cause a vaporisation of coagulation of the tissue, depending on the power density. The hereby caused side effects should be evaluated by a specialised physician before the treatment. Special attention must be paid, to avoid the unintentional coagulation of important structures and tissues during deeper coagulation procedures. The application method should be selected accordingly.

The attending physician has to decide whether or not the given method is appropriate for each individual treatment.

Please refer to chapter *Laser induced fire hazard*.

The selection of the indications suitable for laser treatment in the different special fields is based on their suitability for an endoscopic (and therefore minimally invasive) coagulation or removal of the tissue and is numerous.

Additionally to specific journals and books the collection of single sheets of the ecomed Verlag is recommendable:

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In this manual the basics as well as the different application technique, therapeutic guidelines and the procedures for various special fields, including the appropriate anaesthesia, are described in detail.

Summary of the possible indications for NIR laser according to the therapeutic guidelines (from the source above):

Tip

This summary is given for NIR lasers with a power higher than 25W. The *QuadroStarPRO* may only be used for smaller lesions from the list below.

SCORE:

Indication

- 0 = Contraindication
 1 = not recommended
 2 = possible,
 if other techniques do not work
 3 = possible in competition
 with other techniques
 4 = good indication
 5 = absolute indication,
 another technique is not recommended

procedure

- 0 = impossible
 1 = hardly possible
 2 = hardly possible,
 applicable when no other alternatives possible
 3 = possible
 4 = easily feasible
 5 = very feasible

Technique	Non contact	Contact
<i>Face, oral cavity, pharynx</i>		
Carcinoma (Lip, tongue)	2 / 4	
Adenoids	3 / 4	
Laser assisted uvula plastic (LAUP)	3 / 4	
Tonsillotomy	3 / 4	
Tonsillectomy	2 / 3	
<i>Upper respiratory tract</i>		
Nasal polyp relapse	2 / 3	4 / 5
Nasal concha cauterisation	2 / 4	4 / 5
Nosebleed	5 / 5	
Papillomatosis	4 / 4	4 / 5
<i>Lower respiratory tract</i>		
Tracheal and bronchial fistula (max. 2mm)		4 / 4
Benign tracheal and bronchial stenosis		
- congenital		4 / 5
- scarred		4 / 5
- congenital vascular disease (CVD)	5 / 5	2 / 5
- Papilloma		3 / 4
- Granuloma		3 / 4
Malignant tracheal und bronchial stenosis	4 / 4*	4 / 4*
<i>Thoracic wall and chest cavity</i>		
Thoracic wall		
- Mastectomy		2 / 3
- Thorax wall tumours	3 / 4	3 / 4
- Metastasis	4 / 4	4 / 4 (ITT)
Chest cavity		
- Fistula	3 / 5	
- Decortications		3 / 5
Thoracoscopy		
- Decortications		4 / 5
- Pleurodesis (pneumothorax relapse)	4 / 4	
- Sympathectomy	5 / 4	4 / 5
- Wedge resection	3 / 4	4 / 5

* partly alternating

Gastrointestinal tract		
Esophagotracheal fistula (diameter max. 2mm)		4 / 3
Esophagostenosis		
- Tumorous	5 / 4*	5 / 4*
- Chemical burn		4 / 4
- Barrett's Ulcer	4 / 4	
- congenital		4 / 4
Bleeding		
- Angiodysplasia	5 / 4	
- Ulcers	2 / 3	
Polyp		
- Polyposis coli		3 / 5
- Villous adenoma	3 / 5	
Colorectal carcinoma, pall. Tumour reduction	5 / 4*	5 / 4*
Recanalisation	5 / 4	5 / 4
Abdominal cavity		
Cholesystectomy	3 / 5*	3 / 5*
Appendectomy	3 / 5*	3 / 5*
Adhesiolysis		
- acquired		4 / 5
- inflammatory		4 / 5
Tumour destruction	4 / 5	
- LITT		5 / 5
Lymphadenectomy		4 / 5
Orchidolysis	4 / 4*	4 / 4*
Vagotomy		3 / 4
Abdominal cavity (open)		
Tumours		
- Resection	3 / 4	3 / 4
- Preparation	3 / 4	3 / 4
Parenchymatous Organs		
- Resection	4 / 4	
- Ruptures	4 / 4	3 / 4
- Destruction	4 / 5	
- LITT		5 / 5
Adhesiolysis		
- Congenital		4 / 5
- inflammatory		4 / 5
Lymphadenectomy		4 / 5
Genitourinary system (Gynaecology)		
Condylomata acuminata		
- cornificated epithelium	3 / 5	
- not cornificated epithelium (Mucosa)	5 / 5	4 / 5
Epithelial dysplasia: M. Bowen	2 / 5	
Hysteroscopy		
- benign Endometrial polyp		3 / 5
- submucous Fibroma		4 / 4
- uterine septa		4 / 5
Laparoscopy		
- Endometriosis	4 / 4*	4 / 4*
- Peripubertal adhesion	4 / 5*	4 / 5*
- Polycystic ovaries	4 / 5*	4 / 5*
- Ectopic pregnancy		4 / 5

* partly alternating

Genitourinary tract (gynaecology)		
Condylomata acuminata		
- cornificated epithelium	3 / 5	
- not cornificated epithelium (Mucosa)	5 / 5	4 / 5
Epithelial dysplasia		
- M. Querat	4 / 4	
- M. Bowen	2 / 5	
Bladder tumours		
- Carcinoma in situ	3 / 5	
- Invasive Carcinoma (T ₁)	3 / 5	
Ureter- /urethra stenosis		
- Congenital	4 / 5	
- Acquired		
• benign	3 / 5	
• malignant	3 / 5	
Prostate		
- BPH: transurethral, interstitial		4 / 5 (ITT)
- Prostate carcinoma	2 / 4	2 / 4
- Laparoscopic. Lymphadenectomy		4 / 5
Penis carcinoma	3 / 5	3 / 4
Heminephrectomy	4 / 4	4 / 4
Proctology		
Condylomata acuminata		
- cornificated epithelium	3 / 5	
- not cornificated epithelium (Mucosa)	5 / 5	4 / 5
Epithelial dysplasia: M. Bowen	2 / 5	
Anal fissure		
- Excision		4 / 5
- Coagulation	4 / 5	
Anal marisca		4 / 5
Haemorrhoidectomy		4 / 5
Fistula		
- Congenital		3 / 5**
- Inflammable		3 / 5**
- Tumorous		4 / 5**
Anal stenosis		4 / 5
Pilonidal sinus		4 / 4
Polyp		
- Polyposis coli		3 / 5
- Villous Adenoma	3 / 5	
Skin and skin appendage		
Condylomata acuminata		
- cornificated epithelium	3 / 5	
- not cornificated epithelium (Mucosa)	5 / 5	4 / 5
Epithelial dysplasia: M. Bowen	2 / 5	
Molluscum contagiosum	3 / 4	
Verrucae		
- juvenile	4 / 5	
- palmaris, ungual	3 / 4	
- plantaris	3 / 4	
Benign skin tumours	2 / 4	
cutaneous / Subcutaneous Metastasis (inoperable)	3 / 5	
Basalioma	2 / 5	

*partly alternating

**depending on the extension

Vascular system		
Perforating varixes		4 / 4 (ITT)
Subcutaneous Varixes (AV-Shunts)		4 / 4 (ITT)
congenital: vascular malformation		
- tuberosa	5 / 5***	
- venous, lymphatic, arterial, mixed	5 / 5	5 / 5 (ITT)
M. Osler	5 / 5	
Haemangioma		
- prodromal	3 / 4	
- cutaneous	5 / 4***	
- subcutaneous / organ infestation	5 / 5***	4 / 5 (ITT)
Oncology		
Palliative		
- Metastasis	4 / 4	
- Local recurrences		4 / 4 (ITT)
- Recanalisation	4 / 4*	5 / 4*

* partly alternating

** depending on the extension

***with ice cooling

The **contraindications** have to be considered carefully, according to the various indications in the special fields. Due to the gentler procedure and hemostasis, a laser treatment of the given indications, compared to other methods, is in general more suitable especially for older and multi-morbid patients.

The possible **side effects** that can occur after treatment of the skin or near the skin surface are given in the chapter above *Possible side effects* for the treatment with yellow or green wavelengths (e.g. erythema, oedema, blisters, scars, pigmentary changes). Yet the NIR laser possesses different absorption characteristics and therefore pigmentary changes occur less often (especially the treatment for patients with darker skin type is no contraindication) although occurring scars might be more severe (deeper penetration of the radiation).


These side effects may also occur during or after other treatments, yet they might be undisturbing or even intended.

Accessories

The use of accessories not tested and released by Asclepion Laser Technologies GmbH for the use with this laser device can cause severe damage to the eye or electrical shocks.

Only use accessories and cables provided by Asclepion Laser Technologies for this device.

Important information

We strongly advise against using accessories of other manufacturers. An official testing authority may have certified that a specific accessory unit could be used safely, yet Asclepion Laser Technologies GmbH cannot assume any liability for these products.

On request, Asclepion Laser Technologies will readily send you an updated list of the available accessories.

8 Cleaning, disinfection and sterilisation

8.1 General rules

Device

The laser is operated with line voltage. There is the risk of an electrical shock!

Detach line power cable from the mains socket in order to break the power supply of the system before you proceed with work for cleaning.

Only outer system surfaces may be cleaned by the user.

The accessible system components may be cleaned, using a soft, slightly moistened piece of cloth. A wet piece of cloth must not be used for this purpose.

To remove stubborn contamination, a **soft** cleaner solution or disinfectant may be applied. Refrain from the use of aggressive disinfectants or scouring agents.

Focusing Handpiece and Scanner

The contact surface of the handpiece contacts the client's skin during the treatment. There is the risk of microorganisms being transferred.

Handpiece parts must be cleaned and disinfected immediately, with the device switched off, following use and before application on a client.

If it should be necessary to clean the device during a treatment the device must be switched to the STANDBY mode before starting cleaning works! Do not use alcohol or isopropyl alcohol near the point of treatment because this causes risk of fire if the laser is activated afterwards!

CAUTION

Dirt settling on an optical surface may burn in when exposed to laser radiation and lead to a reduced output and an inhomogeneous beam causing side effects such as reddening, possibly blistering, incrustation and seldom scars!

Ensure that all accessible optical surfaces are not damaged and clean after the cleaning process and remain clean during the treatment.

The procedure instructions given in the following sections are to be followed.

Fibres and application parts without optics

The fibres and application parts contact the tissue of the patient during the treatment. There is the risk of transmitting microorganisms. Some of the fibres and application parts are marked as single use products and should be disposed after the use. It is prohibited to use these products on more than one patient!

Some of the fibres and application parts can be used several times. These have to be thoroughly cleaned, disinfected and sterilised immediately after each treatment. Before each treatment the sterile packaging and the instruments have to be inspected for its intactness and usability.

For the processing of resterilizable fibres there is a separate detailed instruction, which is delivered with the fibers and has to be observed carefully.

If it should be necessary to clean the fibre or other components during a treatment the device must be switched to the STANDBY state before starting cleaning works! Do not use alcohol or isopropyl alcohol near the point of treatment because this causes risk of fire if the laser is activated afterwards!

CAUTION


Dirt settling on the end of the fibre may burn in when exposed to laser radiation and thus causing heating and fire hazard when working with a non contact method.

The end of the fibre must not be damaged and must stay clean during the procedure when working with a non contact method. In case of damage or contamination of the end of the fibre it has to be cleaned or, if necessary, revised.

The procedure instructions given in the following sections are to be followed.

8.2 General advice for processing of the application parts

Application parts are: focusing handpieces, scanner, application parts without optics (for reusable fibres see instruction delivered with the fibres)

Important information

You are prohibited from cleaning a handpiece base body, which contains optics or electronics, under running water. This may destroy the internal optics and electronics. These handpieces cannot be sterilized.

Generally, neutral cleaning agents (e.g. Edisonite Super) are to be used! (The parts are damaged by strongly alkaline solutions!)

Special care is required when cleaning and disinfecting long and narrow cannulas and blind holes, if existing. Please observe the special recommendations of the following section, too.

Frequent cleaning and disinfection have low influence on the accessories. The colour may become lighter or changes but this does not affect the use.

The end of the product service life is normally determined by wear, tear and damage caused by its use.

If the handpieces are not in use for a longer period they must be protected against contamination.

INSTRUCTIONS

Workplace

The workplace must be clean and disinfected.

Preparations for cleaning

The cleaning of the accessories should be performed as soon as possible after their use. For this purpose they are to be removed from the device.

Disassemble accessory units, if required, by following the instructions provided in the next chapter.

Cleaning: manual

Equipment: neutral cleaning agent (e.g. Edisonite Super); brush, ultrasound device (if necessary), wet, lint-free cloth/swabs; soft brush; Q-tip; microfiber wipe delivered by Asclepion or similar lens cleaning tissues; alcohol or isopropyl alcohol

1. Procedure for the handpiece without optics, electronics, cannulas and similar mechanical application parts

The accessible surfaces must be thoroughly cleaned. Remove the dirt by cleaning the accessories part under running water (temp. <30°C). Apply a neutral cleaning agent with a brush on all the surfaces. Use a suitable brush for the cleaning of the cannulas or blind holes. Make sure that the whole surface has been cleaned properly. If necessary use an ultrasound cleaning device. Make sure that the cannulas are cleaned with running water and repeatedly fill and empty blind holes with water.

2. Procedure for the body of the handpiece with integrated optics and/ or electronics

The accessible surfaces of the handpiece body are cleaned by using a soft, slightly damp cloth. The cloth must not be wet or even dripping.

You can use neutral cleaning agent for stubborn dirt. When doing this make sure that possibly existing optical surfaces in the body interior are left out in this process!

3. Procedure for optical surfaces:

The accessible optical surfaces must be absolutely clean and lint-free. Clean them only, if this is not the case.

Observe the optical surfaces sideways in order to assess the pollution, as the dirt is easy to see this way.

First try to remove the dirt by blowing it away or by using a brush. If required, the surface is breathed on and then cleaned with a microfiber (as indicated above). The Q-tip and the microfiber wipe above the Q-tip can be used for areas that are difficult to access. Moisten the wipe with a little bit of alcohol and repeat the process for stubborn dirt.

Make sure that these cleaning agents only get in contact with glass or metal surfaces of the handpiece:

The optical surface must not get scratches and afterwards it must be lint-free and free of streaks.

Cleaning: automatic	Only possible for mechanical components without optics or electronics 1. Place the components in a suitable position, allowing the water to drain from cannula openings and blind holes. 2. Use the standard cycle, wash for at least 5 and rinse for at least 3 minutes. 3. When taking out the pieces check for visible dirt, especially in blind holes and cannulas. If necessary EXPERT the cycle or clean the pieces manually.
Disinfection	The accessible outer surfaces are to be disinfected by means of a disinfection cloth (e.g. Kodan cloth). When doing this, follow the manufacturer's instructions.
Sterilization	Only possible for mechanical parts without optics or electronics Vacuum autoclave, at least 5 minutes at 134°C. The temperature must not be higher than 150°C.
Package for sterilized parts	The packaging must have an appropriate size, so that the fastener is not under pressure. No other specific requirements.
Storage and transport	The parts must be stored in a clean condition and protected against contamination. Ensure that the parameters of the correspondingly admissible ambient conditions indicated in chapter 3 <i>Technical Data</i> are given. The according storage conditions for sterile packages must be met.

The manufacturer has validated that the above instructions are **SUITED** for preparing the accessory parts for their reuse.

The contact surface of the handpiece contacts the client's skin during the treatment.
There is the risk of microorganisms being transferred.
The person who cares about cleaning has all the responsibility that the actually performed cleaning and disinfecting works by means of the equipment and materials and with the staff involved achieve the desired results.

To guarantee this, validations and routine monitoring of the procedures are normally required. The effect and adverse consequences of each deviation from the instructions given should also be evaluated.

8.3 Special recommendations for cleaning and disinfecting the handpieces and scanners

The outer surface of the focusing handpiece and scanner can be cleaned using a slightly moistened cloth as described above.

Rubber bands have to be exchanged if they are worn out.

The **handpieces without optics** have autoclavable cannulas and provides better grip and guidance of the fibres only. They have no electronics as well and can be sterilized also if necessary.


Figure 19: Handpiece without optics with different cannulas

All components of the handpiece are easily removed by pulling and put together. No screwing is necessary (except for the transfer fibre of the device, which usually does not have to be removed).

All glass surfaces (spot size optics, optics adapter and end of fibre in the SMA plug) have also to be kept clean.

When using the **focusing handpiece** the **spacer** contacts the – usually intact – skin of the patient and therefore has to be cleaned and disinfected after each treatment. The spacer has to be detached from the handpiece in order to be cleaned and disinfected. As it is a metal component without any optics or electronics autoclaving is possible, if necessary.

When the spacer is removed from the handpiece in order to be cleaned or disinfected, the accessible surface of the lens inside the spot size optics can be checked for dirt. If this surface is polluted it has to be cleaned using a Q-tip as described above.

The **SMA connector** (of the transfer fibre) should not be detached from the device after the installation, as a change of the spot size and the handpieces can be achieved by changing the optics adapter. For this reason the end of the fibre will most likely not be polluted. In case of a changing of the fibre, special attention has to be paid to ensure a clean fibre end before reattaching the SMA connector.

As long as the same spot size is used and the **optics adapter** is not detached from the body or the handpiece it can be considered as being clean. If the optics adapter has to be removed, the optical surface of the lens inside has to be checked for dirt and if necessary cleaned before reattaching it to the device.


Figure 20: handpiece with transfer fibre

When using the **scanner**, the cleaning procedures and requirements for the SMA connector and the optics adapter are the same as when using the handpiece (as described above).

The cooled **sapphire window and the surrounding metal surface** are in contact with the – usually intact – skin of the patient. Therefore they have to be cleaned and disinfected after each treatment. During the treatment the sapphire window may fog up and the small water drops can reduce the fluence of the laser beam. This can lead to a less effective or even ineffective treatment. Therefore the upper side of the sapphire window has to be cleaned as required with the microfiber wipe provided by Asclepion Laser Technologies or a similar tissue. If the optical surface cannot be cleaned with a wipe, a swab moistened with alcohol or isopropyl alcohol can be used.

The glass surface of the scanner optics usually remains clean, yet it has to be checked for dirt or pollutions after each treatment and if necessary cleaned with a microfiber wipe or a similar tissue.


Figure 21: Scanner with integrated skin cooling

8.4 Special recommendations for cleaning, disinfecting and sterilising the bare fibre and other special application sets

For the indicated re-usable, re-sterilizable fibres the delivered instruction are to be observed, which describe the processing and usage in detail.

Special application sets (EndoSet, Lipo-Set) are intended – in accordance with their labeling – for single use only. The sterility after processing cannot be guaranteed.

The same will apply for other fibers and application parts, which are labeled for single use.

9 Fault messages and fault removal

The laser is operated with line voltage and emits radiation in the visible and the infrared range. There is the risk of an electrical shock and damage to the eye!
You are strongly advised against trying to do repair works at the device yourself that go beyond the instructions given on the display!
Before the change of fuses (at the rear side of the device, below the main switch) the device has to be disconnected from the mains by pulling out the power plug.

9.1 General information

If the device does not respond to turning the key switch, make sure that the following requirements are met:

1. The device is connected to the line power supply and fuses of the room / house are in good order
2. The power switch of the device is switched on (position "I")
3. Fuses of the device are in good order (before device has to be disconnected from the mains by pulling out the power plug)

In case of a general fault or operating error, the display will show a fault message or a note how to solve the problem (which has to be confirmed by the button at the top left, if at the top left of the message OK is displayed). E.g. the one shown below:


Figure 22: Display fault message (example)

In the case of safety-relevant faults, the device will shut down within milliseconds. A display of a fault message is impossible. If this happens or any other fault appears that has not been described in this manual, disconnect the device from line power by unplugging the power plug from the power outlet.

This chapter provides information on possible causes of any malfunction of the QuadroStarPRO system. Faults can be removed by the operating personnel except of those specially referred to as absolutely requiring repair by technicians employed with or authorized by Asclepion Laser Technologies GmbH.

Improper use or maintenance of the laser system may void the warranty granted by Asclepion Laser Technologies GmbH. Before starting any repair or if you have questions that have not been answered by this manual, contact your service partner (see chapter *Customer service*).

9.2 Fault messages and fault removal

The *QuadroStarPRO* is provided with alarm systems. If they are triggered, the system will switch to the safe mode and stop the emission of laser radiation. Please try to remove the fault yourself by following the instructions given in the table below. If the error persists, contact Asclepion Laser Technologies GmbH or a service partner expressly authorized by Asclepion Laser Technologies GmbH.

Before you contact the Technical Service of Asclepion Laser Technologies GmbH, **please make a note of the fault message and the error code displayed.**

Fault message	Explanation / Corrective action
Communication error (all 50 XXX errors)	The communication is interrupted. Restart the device. Operation can be continued. Please contact the Technical Service if the error occurs again.
Temperature error (all 150 XXX errors)	The operating temperature is too high. Check if the device is sufficiently ventilated. Leave the device to cool in STANDBY. Restart the device after some time. Operation can be continued. Please contact the Technical Service if the error occurs again. The operating temperature is too low. Leave the device to acclimate until the minimal working temperature is reached. Operation can be continued. Please contact the Technical Service if the error occurs again.
Foot switch (400 XXX errors)	Please release the foot switch (do not activate the foot switch too early, wait until the display shows the yellow warning sign) !
Door interlock (600 601 errors)	The door interlock contacts are open or the door interlock plug is not connected. If the door has been connected to the door interlock, door is open. Check the door interlock plug. If the door has been connected to the door interlock, close the door.
Handpiece error (all 110X errors)	Electrical fibre plug of the transfer fibre not connected or damaged. Check the handpiece connection. Restart the device. Please contact the Technical Service if the error occurs again.
Unknown applicator (XXX error)	Electrical fibre plug not connected or damaged. Check the fibre connection. Restart the device. Please contact the Technical Service if the error occurs again.
Fibre not connected (XXX error)	SMA - fibre plug not connected or not correctly connected or damaged. Check the SMA - fibre plug. The fibre plug has to be fitted to the stop into the SMA-socket of the device and the screwcap has to be screwed clockwise firmly. Restart the device. Please contact the Technical Service if the error occurs again.
Scanner (XXX error)	Malfunction of the scanner or SCAN mode. Remove the cable of the scanner and connect it once more. Restart the device. Please contact the Technical Service if the error occurs again.

Possible malfunctions

Malfunction	Explanation / corrective action
No power supply	Verify that the power cable has been connected to a power outlet. Is the power switch set to position 'I'? Check both fuses of the device (pull the mains plug before !) Check the line fuses of the room/ domestic power supply system.
No pilot beam	Device is not in READY. Pilot is switched OFF in <i>TOOLS</i> or at too low intensity. Distance of the focusing handpiece too big (pilot beam is weaker if the spot is large). Application parts are not connected correctly. Check that the SMA-plug of the fibre is connected correctly. Optics of the application part is not in good order. Broken fibre (to see at Bare Fibres only: red light along the fibre)
Function keys of the display do not react (or only sporadically)	Faulty calibration of the touch screen The calibration is to be performed by an authorized service engineer. Therefore, inform your distributor or the Technical Service of Asclepion Laser Technologies GmbH (see chapter <i>Customer service</i>).
No laser emission with the foot switch pressed	Verify that the connecting cable of the foot switch is properly connected. The device is in <i>STANDBY</i>. Press the <i>READY</i> key to switch to the <i>READY</i> mode.

10 Customer service

The laser is operated with line voltage and emits radiation in the visible and the infrared range: There is the risk of an electrical shock and damage to the eye!

Never open the device and attempt to repair the device yourself if the laser does not work properly. Only service personnel from or trained and authorised by Asclepion Laser Technologies are allowed to open the device and/ or to perform service work. The same applies to the mandatory annual safety checks, including the calibration.

The operation, maintenance, test and calibration instructions given in this user manual are to be followed.


If you need assistance by the Technical Service, please contact your distributor or us via the address or phone number given below. When contacting your service partner, always keep the model and the serial number of the device at hand. These numbers are specified on the type plate located at the rear panel of the laser device.

Asclepion Laser Technologies GmbH
Service
Brüsseler Str. 10
D - 07747 JENA
Germany
Tel.: +49 (0) 3641 / 7700 - 401
Fax.: +49 (0) 3641 / 7700 - 402
E-Mail: service@asclepion.com

Please note that the manufacturer, installer and importer will not consider themselves responsible for the effects of servicing on safety, reliability and performance of the device unless the following requirements are met:

- installation, extensions, readjustments, modifications or repairs have been performed by authorized persons,
- the electric installation conditions of the laser room meet the requirements of IEC 64 and other national or local standards, as amended, and
- the device has been used in accordance with the instructions given in this manual.

11 Regular maintenance, safety checks and calibration

Important Information

It is absolutely necessary that a preventive maintenance including a safety control of the device is performed **annually** by service personnel from or trained and authorised by Asclepion Laser Technologies in order to operate the device in an optimal condition. This annual check is a requirement for the granting of an extended warranty (after you have registered the device on our homepage).


CAUTION

A too high output power of the laser can cause stronger side effects.

The output power of the therapeutic laser must be checked by a service technician of Asclepion Laser Technologies GmbH or a qualified service technician of the distributor within the technical safety check at least **once a year**. If required the device is also to be calibrated.

The results of the annual safety check have to be documented in the medical devices book.

After a successful calibration the following (test) treatments should be performed with special caution as the fluence value for the identical adjusted fluence parameter might have increased (to the original condition). The service personnel will tell you if and how much the calibration had to be readjusted.

Important Information

The regularly and carefully performed cleaning and maintenance as described in chapter *Cleaning, Disinfection and Sterilisation* and in this chapter ensures a long useful life if the device. Yet the device also contains components that are subject to wear and tear and are not included in the warranty. A list of the wear parts of this device can be requested from the customer service.

Defective or malfunctioning components have to be replaced immediately!

Use only spare parts provided by Asclepion Laser Technologies. When using components of external manufacturers the compliance of the safety conditions cannot be guaranteed, even when the components seem to be equivalent to the ones from Asclepion Laser technologies.

Any modifications of the device are strictly prohibited.

11.1 Routine maintenance

The routine maintenance procedures described below can be performed by the personnel of the clinic or hospital.

General

Inspection of the outer surfaces of the device

Check if all cables are intact and connected firmly.

Frequency: weekly

To be checked by: personnel of clinic or hospital

Check of the cooling system

Check if all ventilation grids of the cooling system are free.

Frequency: weekly

To be checked by: personnel of clinic or hospital

Check of the handpiece and other application parts

Check if the application parts are in perfect state.

Frequency: before each use

To be checked by: operator of the device

11.2 Regular safety check

SCOPE OF ANNUAL SAFETY CHECK

We strongly recommend performing the safety check annually and including at least the tests mentioned below. The tests should be run in the specified order:

- visual inspection of laser device and accessories
- measurement of protective ground wire resistance
- earth leakage test under normal conditions
- measurement of the actual output power of the laser modules at the handpiece or other application parts; calibration of laser power, if necessary
- functional test.

Tip

To inform yourself about the extent and performance of the safety checks, please contact your distributor or directly Asclepion Laser Technologies GmbH.

Additional Information


Upon request, the manufacturer will provide circuit diagrams, component parts lists, descriptions, calibration instructions, or other information not yet contained in this User Manual to assist appropriately qualified technical staff to repair those parts of the device that have been designated by the manufacturer as repairable. In this context, "appropriately qualified technical staff" means staff that attended a manufacturer's service training course on this device and that was authorized to repair it.

12 Disposal

The device is to be disposed of according to regulation 2002/96/EC of the European Council for Waste Electric and Electronic Equipment [WEEE] or the specific national regulations on waste electric and electronic devices, as amended.

Please consult our Technical Service if you have any questions regarding disposal.

13 Manufacturer's declaration of EMC compliance

Important Information

The use of accessory parts, transformers or wiring details other than specified hereafter may result in increased spurious emission and /or reduced noise immunity of the device.

Accessory part/ designation	Item Number	Length/ phys. dimensions
Footswitch with cable	660806008	< 3,0 m
cold appliance connection cable	550704001	< 3,0 m
Scanner	4048	< 3,0 m

Important Information

The device must not be stacked / arranged directly besides or with other devices. If used directly besides or with other devices the device has to be observed to control his intended use in the arrangement in question.

3.1 Emitted electromagnetic interference

The *QuadroStarPRO* is solely designed for operation in an electromagnetic environment as specified below. The Customer and/or Owner/Operator should ensure that the *QuadroStarPRO* will operate in an electromagnetic environment according to specifications quoted hereafter.

Measurement of the emitted interference	Compliance	Guideline features regarding electromagnetic environment
High-frequency spurious emission under CISPR 11	Group 1	The <i>QuadroStarPRO</i> uses HF energy solely for its own function. For this reason, spurious emission of high-frequency interferences is very small and other electronic equipment located in the surrounding area is unlikely to be disturbed.
High-frequency spurious emission under CISPR 11	Class B	The <i>QuadroStarPRO</i> is designed for use in all institutions, including residential quarters and any operating site with a direct connection to the public power network that supplies buildings for residential purposes.
Emission of spurious harmonics under IEC 61000-3-2	Class A	For harmonics: The <i>QuadroStarPRO</i> is designed for use in all institutions, including residential quarters and any operating site with a direct connection to the public power network that supplies buildings for residential purposes.
Spurious emission of voltage variations/flicker under IEC 61000-3-3	Complies	For flicker: The <i>QuadroStarPRO</i> is designed for use in all institutions, including residential quarters and any operating site with a direct connection to the public power network that supplies buildings for residential purposes.

13.2 Electromagnetic Immunity

The *QuadroStarPRO* is intended for operation in a specified electromagnetic environment. It is for the Customer and/or Owner/Operator to ensure that the *QuadroStarPRO* operates in an electromagnetic environment as described hereafter.

Test for noise immunity	Test level under IEC 60601-1-2	Compliance level	Guideline features regarding electromagnetic environment
Discharge of static electricity (ESD) under DIN EN 61000-4-2	+/- 6 kV contact discharge +/- 8 kV air discharge	+/- 6kV +/- 8kV	Flooring should consist of wood, concrete or a finish coating of ceramic tiles. Where flooring is provided with synthetic material, relative air humidity must not be less than 30%.
Fast transient electrical interferences/bursts under IEC 61000-4-4	+/- 2 kV for line power cables +/- 1 kV for input and output cables	+/- 2kV	The quality of supply voltage should be that of a typical business and/or hospital environment.
Surge voltages (surges) under IEC 61000-4-5	+/- 1 kV voltage outer conductor-outerconductor +/- 2 kV voltage outer conductor-earthing	+/- 1kV +/- 2kV	The quality of supply voltage should be that of a typical business and/or hospital environment.
Voltage dips, short-time breaks and variance in supply voltage under IEC 61000-4-11	< 5% UT for ½-period (> 95% dip) 40% UT for 5 periods (60% dip) 70% UT for 25 periods (30% dip)	< 5 % 40% 70%	The quality of supply voltage should be that of a typical business and/or hospital environment. Users requiring uninterrupted function of the system, including in the case of a break in energy supply, are recommended to use an uninterruptible power supply module or a battery backup for power supply.
	< 5% UT for 5 s (> 95% dip)	< 5 %	
Magnetic field for supply frequency (50/60Hz) under IEC 61000-4-8	3 A/m	3 A/m	Magnetic field strength should correspond to the typical value of business and hospital environments.
Additional remark:	U _T designates the line input AC voltage prior to application of test levels * - untested, because cabling length is less than 3 m		

The *QuadroStarPRO* is designed for operation in an electromagnetic environment that conforms to specifications. It is for the Customer and/or Owner/Operator to ensure that the *QuadroStarPRO* operates in an electromagnetic environment as described below.

Test of noise immunity	Test level under IEC 60601-1-2	Compliance level	Guideline features regarding electromagnetic environment
Conducted HF interferences under IEC 61000-4-3	3 Veff 150 kHz to 80 MHz	3V	Portable and non-stationary radio equipment should not be operated at a distance from the system, including related cabling, that is smaller than the recommended safety distance. This value is calculated according to an equation formula that is adapted for the frequency of emission. Recommended safety distance: $d = \left[\frac{3.5}{3V} \right] \sqrt{P}$
Emitted HF interferences under IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$ for 80 MHz to 800 MHz $d = \left[\frac{7}{3} \right] \sqrt{P}$ for 800 MHz to 2.5 GHz where P designates the nominal emitter power in watts (W) according to specifications of the emitter manufacturer and d stands for the recommended safety distance in meters (m). According to on-site test results ^a , the field strength of stationary radio emitters should be below the level of compliance at all frequencies. ^b Interferences may be expected in the proximity of equipment with a pictogram as follows. 

^a The field strength of stationary emitters such as a base station of a cell phone (cellular/cordless), mobile land radio devices, amateur radio equipment, AM and FM broadcasting and television broadcasting cannot be pre-determined with theoretical precision. To analyze a given electromagnetic environment for stationary transmitters, a special site study should be considered. Where readings measured for field strength at the point of use are found to exceed the above quoted levels of compliance, the *QuadroStarPRO* should be kept under observation, in order to ensure intended use and conforming operation. Where fitness for use is found to be non-conforming, additional measures may be required, e.g. a modified alignment position or relocation to another site.

^b For the frequency range from 150 kHz to 80 MHz, field strength should be less than 3 V/m.

Additional note: These guidelines may not be applicable in all cases. Absorption and reflection caused by buildings, objects and people may have an influence on the propagation of electromagnetic quantities.

13.3 Recommended Safety Distances

Recommended safety distance values between portable or mobile RF telecommunication devices on the one and the *QuadroStarPRO* on the other hand (under DINEN60601-1-2):

Transmission frequency	150 kHz to 80 MHz	150 kHz to 800 MHz	800 MHz to 2.5 GHz
Equation	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
Nominal power of transmitter (watts)	Safety distance (in meters)	Safety distance (in meters)	Safety distance (in meters)
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,70	3,70	7,37
100	11,70	11,70	23,30

For transmitters with a maximum nominal power not contained in the table above, the distance value can be determined using the equation formula that relates to the particular column, where P designates the maximum nominal power of the transmitter in watts (W) according to specifications by the transmitter manufacturer.

Additional note:

These guidelines may not be applicable in all cases. Absorption and reflection caused by structures, objects and people may have an influence on the propagation of electromagnetic quantities.

14 EC declaration of conformity as medical device - Copy

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
Jena, May 21st 2013

Asclepion Laser Technologies GmbH
i.V.

Pietro Amatulli
Technical Director

Directive 93/42/EEC, all modifications included

15 Appendix: Cosmetic treatment with the infrared QuadroStarPRO

The device confirms with the requirements of a medical device (see chapter *EC declaration of conformity as medical device*) solely for the intended use given in chapter 1.3.
The cosmetic application stated in the next chapter does not apply the evaluation of conformity according to the guideline 93/42/EEC of the certificate of chapter 14.

Moreover the device can also be used for cosmetic applications by a doctor or under responsibility by a doctor.

For this cosmetic application all instructions of this manual apply including warnings for the safe use, operation and maintenance with following supplements:

15.1 Intended use

The relatively unspecific warming by the laser radiation of the infrared version of the QuadroStarPRO can be used for the following cosmetic application also – preferably with the wavelength of 980nm:

- temporary increase of clear nail in clients with onychomycosis (e.g. Trichophyton rubrum and T. menisporum, and/or Candida albicans, etc.).

WARNING:
The incorrect using of this device can cause serious danger to health. Therefore, all warnings, instructions and information of this User Manual are to be observed under all circumstances!

Because this is a cosmetic application the clients have to be healthy besides the onychomycosis. Further restrictions have to be taken from chapter 1.3.2. If there is the slightest doubt if a contraindication applies or not, refuse to do the application until there is the diagnosis by a doctor, that there is no contraindication.

Additional contraindication: clients who have a dysfunction of pain feeling have to be excluded

15.2 Changes in other chapters


The following supplements or ascertainment apply to the cosmetic treatment of nails:

Spezielles Anwendungsteil:	Focusing handpiece for treatment of nails with a spot size of 3.5mm. No differences to the other focusing handpiece described above relating to installation, cleaning/ disinfection.
----------------------------	--

Application:	All basic principles described in chapter Application can be applied. For this treatment of nails the <u>infrared</u> version of the QuadroStarPRO with <u>focusing handpiece</u> is used.
--------------	---

For application the following technique is used:

Application technique: Using of the spot-by-spot technique described already for skin treatments, adjusted to the nail treatment:

The spots are set in serpentine pattern with some overlap back and forth, than top down. About 2mm around the nail is included in the treatment.

Method 1: in constant moving with high repetition rate and low parameters scan the nail until a clear warming is felt, but no pain! If pain comes up STOP immediately and cool! The nail bed must not be damaged! Duration dependent on nail size about 1.5 to 2 minutes.

Method 2: with low repetition rate to set spot by spot next each other. About 3 passes until clear warming is felt, but no pain! If pain comes up STOP immediately and cool! The nail bed must not be damaged! Duration dependent on nail size about 1.5 to 2 minutes.

Follow-up treatments: multiple treatments are necessary with an interval of a few weeks. Longer interval, if the nail grows slowly. Recurrences are possible as with other treatments. The clients have to be informed. Than again a new treatment is possible.

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Instructions for reprocessing of focusing handpieces of the laser QuadroStar

WARNING

The laser device emits strong radiation. That may cause serious damage if the device and the handpieces are used and/ or reprocessed improperly. Moreover there is the danger of microorganism being transferred.

Therefore, always observe the user manual of the laser, these instructions for reprocessing, the requirements of hygiene, if medical products are reprocessed, in your country and your facility as well as the user manuals of the manufacturers of cleaning/ disinfection products and devices as well as of the devices for sterilization! The reprocessing shall only be done by substances and devices, which are CE approved and by staff qualified and authorized to do so.

Valid for: Set focusing handpiece (labeled at spot size optics: **577/532 Ø 1mm***)

optional available spot size optics 0.5mm

optional available spot size optics 1.5mm

optional available spot size optics 2.8mm

*labeling is for the optional spot size optics 0.5/ 1.5 or 2.8mm correspondingly

Set focusing handpiece (labeled at spot size optics: **940/980 Ø 3.5mm***)

item no.: 4042

item no.: 4043


item no.: 4045

item no.: 4046

item no.: 4057

Limitation of reprocessing: Frequent reprocessing has minor effect on the handpiece. The color may become lighter or changes but this does not affect the use. The end of the product service life is normally determined by wear, tear and damage caused by its use.

The handpiece consists of the optics adapter (connected to the transfer fiber), the spot size optics and the distance holder. Only the distance holder can be sterilized, if necessary (metal part, there is no optics and no electronics).

The distance holder contacts the - uninjured - skin of the patient.

The handpiece is delivered not disinfected and the distance holder not sterilized. Before first use and immediately after each treatment the handpiece has to be clean and disinfected, if necessary the distance holder has to be sterile. To do this, proceed as follows:

Instructions for reprocessing	
Work place	The workplace must be clean and disinfected. Two areas have to be prepared: one, which is be used to put down the dirty handpiece, and one to put down the clean and disinfected handpiece. The reprocessing has to be performed in the rooms/ areas intended for these work only.
Preparations for reprocessing	The reprocessing should be performed before first use and within 30 minutes after each use. Use disposable gloves; they have to be liquid-tight. First the device has to be switched off by key switch and main switch, Than spot size optics has to be removed from the optics adapter and the optics adapter has to be put into the handpiece holder of the device <u>carefully</u>. The spot size optics has to be removed from the distance holder and put down at the "dirty" area of the work place. The distance holder has to be put into a container with cold water. Equipment: moistened, lint-free cloth; Kodan® or Bacilloi® wipes; surgical brush; small surgical bottlebrush (no metal brush!); neutral cleaning agent Edisonite® Super; ultrasound cleaning device; Q-tips (lint-free, <u>non-impregnated</u> cotton bud); water with quality of drinking water alcohol or isopropyl alcohol
Instructions for the SPOT SIZE OPTICS	
Cleaning: manual	Machine-based cleaning or cleaning under running water is not allowed for the spot size optics! The outer surface is cleaned using a soft, <u>slightly damp</u> cloth. The cloth must not be wet or even dripping! For stubborn dirt Kodan® or Bacilloi® wipes are used.

Instructions for reprocessing of focusing handpieces of the laser QuadroStar

Additional instructions for the SPOT SIZE OPTICS	
Cleaning: manual	The surface inside keeps clean normally and has to be cleaned only, if this is not the case. To do this, a <u>slightly damp</u> Q tip is used. Make sure that optical surfaces in the interior are left out in this process! Procedure for optical surfaces: The accessible optical surfaces must be absolutely clean and lint-free. After cleaning of the other surfaces, the optical surfaces (internal) are checked for cleanliness by looking through against a window or another light source. If necessary, first try to remove the dirt using a dry Q-tip. Doing this clean from central to outwards by spirals. In case of stubborn dirt moisten the Q-tip with a <u>little bit</u> of alcohol or isopropyl (should not be wet!) and clean from central to outwards by spirals. Under no circumstances do not use these inflammatory substances before the laser is activated or during laser treatment - because this causes a risk of fire and explosion! The optical surface must not get scratches and afterwards it must be lint-free and clean. Dirt settling on the optics may burn in when exposed to laser radiation. This can lead to mistakes and side effects! Check the parts for remaining dirt after drying (if necessary repeat cleaning) and damages. <u>Only impeccable and clean parts can be used further.</u>
Disinfektion	<u>Machine-based disinfection is not allowed for the spot size optics!</u> The outer surface is disinfected by Kodan® or Bacillo® wipes.
Instructions for the DISTANCE HOLDER:	
Cleaning: manual Also for pre-cleaning for machine-based cleaning	The accessible surfaces have to be cleaned carefully. Remove the dirt from distance holder by cleaning under running water (temp. lower than 30°C / 86°F; 3 min. under running water). Clean all surfaces with Edisonite® Super and a brush. Use a suitable bottle brush for the inner surfaces. <u>Only if a machine-based cleaning is not available</u>, clean afterwards with ultrasound device and Edisonite® Super. The cleaning agent must be rinsed thoroughly and the part has to be let dry. Check the parts for remaining dirt (if necessary repeat cleaning) and damages. <u>Only impeccable and clean parts can be used further.</u>
After machine-based cleaning (recommended)	1. Place the distance holder in a suitable position, allowing the water to drain. 2. Use the standard cycle, wash for at least five minutes and rinse for at least 3 minutes. 3. When taking out the parts after drying check for damages and visible dirt – repeat cleaning cycle if necessary. <u>Only impeccable and clean parts can be used further.</u>
Disinfektion	After automated cleaning always a thermal disinfection is used (temperature 93°C / 200°F). After manual cleaning the distance holder is to be disinfected by means of Kodan® oder Bacillo® wipes.
Sterilization	<u>Only distance holders can be sterilized.</u> If necessary, the sterilization is realized using humid heat after disinfection: in a high pressure steam sterilizer at a temperature of 134°C for at least 5 min. or in a steam sterilizer at a temperature of 121°C for at least 20 min.
Packing	Sterile packing is not necessary, because the use of the focusing handpiece under sterile conditions is not necessary and not possible. The sterilization which may be necessary is done due to destroy microorganisms, which can attach at the distance holder from the last treatment (e.g. after treatment of warts).

Instructions for reprocessing of focusing handpieces of the laser QuadroStar

Instructions for the SPOT SIZE OPTICS and DISTANCE HOLDER:	
Reassembling	Immediately after drying the handpiece is reassembled by plugging the distance holder at the spot size optics and this at the optics adapter. It is assumed that the optics adapter is clean, which is normally the case. If the optics adapter is not clean (including optics), it has to be cleaned as described for the spot size optics. It has to be observed, that all parts are connected correctly, i.e. there are no gaps between the parts and the connection is tight. <u>Only impeccable and clean handpieces can be used further.</u>
Storage and transportation	If the handpiece shall not be used for a longer time, it has to be stored fully mounted*, clean and dry in the delivered packing box and - if necessary - to transported. The parameters for storage and transport stated in the chapter <i>Technical Data</i> of the user manual of the device have to be observed. *Fully mounted means, that the handpiece is connected with the transfer fiber via optics adapter. The transfer fiber has to be unscrewed from the back of the device (SMA plug and electrical plug for fiber recognition) and the SMA plug has to be protected by the delivered small cap. Before the transfer fiber is connected again to the device, it has to be checked, that the fiber end in the SMA plug is clean; if it is not, a cleaning has to be done as described for the optics of the spot size optics. Furthermore it has to be checked, that all connectors are tightly screwed/ plugged and all parts are undamaged and clean: fiber end – optics adapter – spot size optics – distance holder.

If there are parts, which are damaged and contaminated, they have to be cleaned and disinfected as described above. If necessary the distance holder has to be sterilized before disposal. This shall guarantee, that a contamination of the parties is excluded.

The person who cares about reprocessing has all the responsibility that the desired results are reached at the basis of the actually performed reprocessing including the chosen tools/ substances, personnel and equipment. She should have a certified qualification to a sterilization assistant.

To guarantee this, validations and routine monitoring of the procedures are normally required. The effects and adverse consequences of each deviation from the instructions given should also be evaluated.

The manufacturer assumes no liability for immediate or consequently damages or injuries in case of mishandling and inappropriate reprocessing.

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Instructions for reprocessing of The scanner of the laser QuadroStarPRO

WARNING

The laser device emits strong radiation. That may cause serious damage if the device and the scanner are used and/ or reprocessed improperly. Moreover there is the danger of microorganism being transferred.

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
Valid for: Scanner for QuadroStarPRO with wavelength of 577 or 532nm item no.: 4048
labeling at the cord of the scanner (see figure below)

Limitation of reprocessing: Frequent reprocessing has minor effect on the scanner. The color may become lighter or changes but this does not affect the use. The end of the product service life is normally determined by wear, tear and damage caused by its use.

The scanner consists of the optics adapter (connected to the transfer fiber), the basic body with integrated optics and electronics and the sapphire plate, which is cooled and serves as cooling of the upper skin layer.

The scanner contains optics and electronics and may certainly not be sterilized or machine-based cleaned or disinfected! The scanner may certainly not be cleaned or disinfected wet or even under running water!

The sapphire plate contacts the - uninjured – skin of the patient.

The scanner is delivered not disinfected. Before first use and immediately after each treatment the scanner has to be clean and disinfected.

To do this, proceed as follows:

Instructions for reprocessing

Work place	The workplace must be clean and disinfected. Two areas have to be prepared: one, which is be used to put down the dirty scanner, and one to put down the clean and disinfected scanner. The reprocessing has to be performed in the areas intended for these work only.
Preparations for reprocessing	The reprocessing should be performed before first use and within 30 minutes after each use. Use disposable gloves; they have to be liquid-tight. First the device has to be switched off. Than the scanner has to be removed from the optics adapter and the optics adapter has to be put into the handpiece holder of the device <u>carefully</u>. The scanner is put down <u>carefully</u> at the "dirty" area of the work place. For this purpose the laser device should be next to the work place, that the plug of the scanner at the rear side of the device must not be disconnected. Equipment: clean, a bit moistened, lint-free cloth; clean microfiber glasses wipe or silicon-free lens paper; Kodan® or Bacillo® wipes; clean Q-tips (lint-free, <u>non-impregnated</u> cotton bud); alcohol or isopropyl alcohol
Cleaning: manual	<u>Machine-based cleaning or cleaning under running water is not allowed!</u> The outer surface of the basic body including the metal surrounding of the sapphire plate is cleaned using the soft, <u>slightly damp</u> cloth. The cloth must not be wet or even dripping! For stubborn dirt Kodan® or Bacillo® wipes are used. The optics (output lens and sapphire plate) has to be excluded.

Instructions for reprocessing of The scanner of the laser QuadroStarPRO

Cleaning: manual	Procedure for accessible optical surfaces (output lens and sapphire plate): The accessible optical surfaces must be absolutely clean and lint-free. After cleaning of the other surfaces, the optical surfaces are checked visually for cleanliness. The output lens normally keeps clean and is cleaned only, if it is necessary. The sapphire plate is always cleaned and disinfected. First try to remove the dirt using the microfiber glasses wipe. Doing this a slight pressure is used only to prevent the sapphire plate is pushed out of the metal mount. In case of stubborn dirt moisten the Q-tip with <u>a little bit</u> of alcohol or isopropyl (should not be wet!) and clean from central to outwards by spirals. Under no circumstances do not use these inflammatory substances before the laser is activated or during laser treatment - because this causes a risk of fire and explosion! The optical surface must not get scratches and afterwards it must be lint-free and clean. Dirt settling on the optics may burn in when exposed to laser radiation. This can lead to mistakes and side effects! Check the parts for remaining dirt after drying (if necessary repeat cleaning) and damages. <u>Only impeccable and clean parts can be used further.</u>
Disinfection	<u>Machine-based disinfection is not allowed for the scanner!</u> The outer surface is disinfected by Kodan® or Bacillol® wipes.
Sterilization	<u>Sterilization is not allowed for the scanner!</u>
Reassembling	Immediately after drying the scanner is reassembled by plugging the scanner at the optics adapter. It is assumed that the optics adapter is clean, which is normally the case. If the optics adapter is not clean (including optics), it has to be cleaned as described for the spot size optics (please see instructions for reprocessing of the focusing handpiece). It has to be observed, that the basic body of the scanner is connected correctly, i.e. there are no gaps between the parts and the connection is tight. The scanner is <u>carefully</u> pulled down in the holder and can be used further. <u>Only impeccable and clean handpieces can be used further.</u>
Packing	Packing is normally not necessary. If the scanner shall not be used for a longer time, the multifunction scanner plug at the rear side of the device has to be disconnected. The scanner has to be stored clean, dry and protected from dirt in the delivered packing box.
Storage and transport	The scanner has to be stored clean, dry and protected from dirt in the delivered packing box and - if necessary - to transported with care. The parameters for storage and transport stated in the chapter Technical Data of the user manual of the device have to be observed.

If there are parts, which are damaged and contaminated, they have to be cleaned and disinfected as described above before it is shipped for repair. This shall guarantee, that a contamination of third parties is excluded. For shipment the delivered packing box has to be used and additionally to package the box in a safe way, so that no damage at transport is possible.

The manufacturer assumes no liability for the transport.

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To guarantee this, validations and routine monitoring of the procedures are normally required. The effort and adverse consequences of each deviation from the instructions given should also be evaluated.

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

Revision index: F

Identification number: MMo/14/02/10-01

Member code/Year/month/day-sequential number

FB 7302 RAaC QuadPro
V03_MMo140210-1.docxMedical Product: **QuadroStarPro**

Item number: 1490 - 1495

Assessed construction level: pre-series

Revision user's Manual: 1490 99 001_Vs2-10.Febr14-
GA_QuPRO-DE☒ MDD 93/42/EEC☒ IEC/ EN 60601-1:2006☒ IEC/ EN 60 601-1-2:2007☒ IEC/ EN 62304:2006☒ IEC/ EN 62366:2008☒ EN 60 601-2-22:2005☒ EN 60 825-1:2007☐☐ before R&D☒ during R&D☐ after R&D☒ (before serial production)☐ after production



Remarks: The risk management of the device was carried out according to EN ISO 14971:2009. At present, an analysis of the implementation of EN ISO 14971:2012 is performed. For this purpose, it is planned to adapt the PA 730 to perform a GAP analysis to the previous made risk management and adapt the risk analysis.

Related documents: Appendix to form / Anlage zum Formblatt FB 7302 Risk Analysis and control Ident.- no. MMo/14/02/10-02

Drawn up:

Vs.	Description of change	Date	Name	Sign
03	risk analysis acc. EN ISO 14971:2009	Feb. 02, 2014	M.Möller	
02	Revision	Oct. 01, 2013	M.Möller	sgd. Möller
01	First edition: Initial version with risk analysis, first risk evaluation and proposal of risk control measures	Nov. 23, 2012	M.Möller	sgd. Möller

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Case of last modification form: new dokument structure with spread sheet

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
File name: FB 7302 RAaC QuPro V03 MMo140210-01.xls

Risikoanalyse / Risk Analysis										Risikobehandlung / Risk Control							
										Sicherheit Safety		Umsetzung der Maßnahmen Counter Measures Implementation			Restrisikobewertung Residual Risk Evaluation		
Nr. No.	Gefährdung Hazard	Med./techn. Auswirk. Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Einschätzung First evaluation	Maßnahme(n) Counter measure(s)	IEC 60601-1:2005 14.6.2 [TRF]	Angewandte Normen Applied Standards	Konstruktiv/Design Schutzmaß/B safeguard Information	Wirksamkeit geprüft? Effectiveness tested?	Umsetzung geprüft? Implementation checked?	Nachweisdokument(e) Record(s)	Verantwortlicher/ Responsible person	A	Restrisikobewertung Evaluation of residual risk	Risiko-Nutzen-Analyse Risk/Benefit Analysis
HA90	Hazards due to insufficient processes	Unwanted side effects (moderate)	critical	Wrong results from analysis due to missing or wrong tools	occasional		In addition to subclause 6.1 of ISO 14971 suitably validated tools and procedures must be selected and identified to implement each risk control measure. These tools and procedures must be appropriate to assure that each risk control measure satisfactorily reduces the identified risk(s).	IEC 60601-1:2005 14.6.2 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			PA-Risikomanagement See FMECA (this document)	IS	inconceivable		
HA100	Hazards due to insufficient processes	Unwanted side effects (moderate)	critical	No / insufficient requirement specification	occasional		For the PEMS and each of its subsystems (e.g. for a PESS) there must be a documented requirement specification: - SRS (Master-Slave)	IEC 60601-1:2005 14.7 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Software Requirement Specification: FB7312 SRS QuadroPRO V01_IS120626-1 Architecture DesignChart PEMS QuadroPro V01 IS120703-03	IS	inconceivable		
HA110	Hazards due to insufficient processes	Unwanted side effects (moderate)	critical	No specification of architecture	occasional		For the PEMS and each of its subsystems, an architecture must be specified that must satisfy the requirement specification.	IEC 60601-1:2005 14.8 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Software Design Specification: FB 73221 SDS GUI Taurus Main V01 NJ120425-01 FB 73221 SDS Keypad Taurus V01 NJ120425-02 FB 73221 SDS GUI Taurus SerialCOM V01 NJ120425-03 FB 73221 SDS GUI Taurus ModDiode V01 NJ120425-04 FB 73221 SDS QuadroPRO V01 FF121023-01 Appendix Architecture Slave V01 FF121023-02	NJ	inconceivable		
HA120	Hazards due to insufficient processes	Unwanted side effects (moderate)	critical	Failure in Firmware / Software and/or omission of required activities and/or no/faulty/incomplete documentation	occasional		The design was decomposed into following subsystems - Master - Slave For every subsystem a design and test specification must create.	IEC 60601-1:2005 14.9 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Software Requirement Specification (see row above) Software Design Specification (see row above)	IS NJ FF	inconceivable		
HA130	Hazards due to insufficient processes	Unwanted side effects (moderate)	critical	No / insufficient planning of the verification	occasional		Verification has to be done for all functions that implement basic safety, essential performance or risk control measures. The verification must be performed according to the verification plan. The results of the verification activities must be documented.	IEC 60601-1:2005 14.10 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			SW Verification Report(s)	NJ	inconceivable		

2.1.2 Line voltage at internal conductive parts (accessible via openings of enclosure)

No. Hazard	Measures to Avoidance Measures to Effect	5	Technical Evaluation Corrective measures	A	Erstbewertung First evaluation	Maßnahmen (Measures)	EC 60601-1:2005 11.6.2 (TRF)	Angewandte Normen Applied Standards	konstruktiv/Design	Schutzmaßn./Safeguard	Wirksamkeit geprüft? Effectiveness tested?	Umsetzung geprüft? Implementation checked?	Nachweisdokumente (Records)	Verantwortlicher/ Responsible person	A	Restriktionsbewertung Evaluation of residual risk	Risiko-Nutzen-Analyse (Risk/Benefit Analysis)
HA240	Line voltage hazards	Injury as a result of electric shock	disastrous	Wrong design of fire device	occasional	There is a liquid storage chamber built into the device. The input is in a closed hood. According to IEC 60601-1, clause 11.6.2, Type of liquid: Glycol Volume: 30ml (15% of 200ml) Duration of spill: 1min After these procedures, the device is to pass the appropriate dielectric strength and leakage current tests n.a., no handling with liquids in normal use	IEC 60601-1:2005 EN 60601-1:2006	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601	LGA	inconceivable		
HA250	Line voltage hazards	Injury as a result of electric shock	disastrous	Wrong design of the device	occasional	Earthing of liquid system	IEC 60601-1:2005 11.6.3 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601	LGA	inconceivable		
HA260	Line voltage hazards	Injury as a result of electric shock	disastrous	Wrong device handling	occasional	Identify hazard due to incorrect connection of accessible connectors. Following interfaces are identified: - USB interface - Service Interface (Sub-D, PS2) - Mains inlet - Doorlock - Food switch interface - Cooling interface - Laser warning lamp - Handpiece - Fiber detector Shape coding of the connectors. No risk of interchanging of the plugs by users, patients and third parties. It is possible to interchange the plugs of doorlock and laser warn lamp but no hazard for users, patients and third parties. Determination of the IP grade: - IP20 Description of the IP-grade in user manual	IEC 60601-1:2005 15.4.1 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601 Datasheet Verification Report(s)	MMo	inconceivable		
HA270	Line voltage hazards	Injury as a result of electric shock	disastrous	Wrong device handling	hardly conceivable	Identify hazard due to incorrect connection of accessible connectors. Following interfaces are identified: - USB interface - Service Interface (Sub-D, PS2) - Mains inlet - Doorlock - Food switch interface - Cooling interface - Laser warning lamp - Handpiece - Fiber detector Shape coding of the connectors. No risk of interchanging of the plugs by users, patients and third parties. It is possible to interchange the plugs of doorlock and laser warn lamp but no hazard for users, patients and third parties. Determination of the IP grade: - IP20 Description of the IP-grade in user manual	IEC 60601-1:2005 15.4.1 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601 Datasheet Verification Report(s)	MMo	inconceivable		
HA280	Line voltage hazards	Injury as a result of electric shock	disastrous	Entry of liquid(s) into the device. Entry of solid materials into the device.	hardly conceivable	Determination of the IP grade: - IP20 Description of the IP-grade in user manual	IEC 60601-1:2005 11.6.5 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	User Manual	ZI	inconceivable		
2.1.3 Short Circuit of insulation barrier																	
HA290	Hazardous voltage is touchable	Electrical shock	critical	Contact with live voltage mains voltage	hardly conceivable	n.a. because no insulators other than wire insulators are contained in the device. - Protective earthing of enclosure with one MOPP - 2xMOPP to SP/SCP - 2xMOPP to applied parts - Housing only be opened with tools	IEC 60601-1:2005 8.8.4.1 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601	LGA	inconceivable		
HA300	Hazardous voltage is touchable (Power supply)	Fire	critical	Contact with live voltage mains voltage	hardly conceivable	n.a. because no insulators other than wire insulators are contained in the device. - Protective earthing of enclosure with one MOPP - 2xMOPP to SP/SCP - 2xMOPP to applied parts - Housing only be opened with tools	IEC 60601-1:2005 8.8.4.1 (TRF) MD 6342/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Test Report of IEC 60601 Insulation diagram: Isolation diagram QuattroPRO V01 MMo121102	MMo	inconceivable		
HA310	Hazardous voltage is touchable	Burns (minor)	critical	Contact with short-circuit current in device and insulation defect Failure of basic insulation	hardly conceivable	- Protective earthing of enclosure with one MOPP - Insulated enclosure - 2xMOPP to SP/SCP - 2xMOPP to applied parts - Mains fuse	MD 6342/EEC	MD 6342/EEC	Ja	Ja	Ja	Ja	Test Report of IEC 60601 Insulation diagram: Isolation diagram QuattroPRO V01 MMo121102	LGA	inconceivable		

File name: FB 7302 RAaC QuPro V03_MMo140210-01.xls

Hazard / Risk Analysis					Risk Management / Risk Control									
Nr.	Hazard	Med./tech. Auswirkung Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Schweregrad Severity	Maßnahme(n) Counter measure(s)	Angewandte Normen Applied Standards	Sicherheits- Safety	Umsetzung der Maßnahmen Counter Measures Implementation	Residualrisikobewertung Residual Risk Evaluation	Residualrisikobewertung Residual Risk Evaluation	Residualrisikobewertung Residual Risk Evaluation	Risiko-Nutzen-Analyse Risk/Benefit Analysis
HA440	Hazardous voltage is touchable	Injury as a result of electric shock	disastrous	(SFC) Accessible conductive part(s) is/are not sufficient separated from applied part(s) of type B not connected with protective earth.	occasional		The spacer of the handpiece and the scanner are applied parts type B. - The applied part are protectively earthed or - separated from metal accessible parts (enclosure) by 1x MOPP	IEC 60601-1:2005 8.5.2.2 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Insulation diagram: Isolationsdiagramm QuadroPRO V01 MM0121102	MM0	Inconceivable	
HA450							No patient leads with conductive material identified.	IEC 60601-1:2005 8.5.2.3 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja				
HA460	Hazardous voltage is touchable	Injury as a result of electric shock	disastrous	Wrong design of foot switch / wrong IP grade	occasional		There is a foot switch: IP grade: IPX6	IEC 60601-1:2005 15.4.7.3 b) [TRF] MDD 93/42/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja Ja	test report foot switch, see LoCC Medizin FS Steute IEC 60601-1,-2-22_TR TÜV + Medizin FS Steute IEC 60601-1,-2-22_TR TÜV_Ergänzung	Mmo (Steute)	Inconceivable	
HA470	Hazardous voltage is touchable	Injury as a result of electric shock	disastrous	Failure of safety relevant components	occasional		Identification of relevant components of device	IEC 60601-1:2005 4.8 [TRF] 4.8 a) [TRF] 4.8 b) [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	LoCC FB 7331 LoCC QuadroPRO V01 MM0121101-01	MM0	Inconceivable	
HA480							No components with high integrity.	IEC 60601-1:2005 4.9 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja				
HA490	Hazardous voltage is touchable	Injury as a result of electric shock	disastrous	(SFC) Not safety relevant design of the device	hardly conceivable		No connectors or conductors are identified to cause potential hazards in this device, if interrupted, unintended moved or accidentally detached.	IEC 60601-1:2005 8.1 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Test Report of IEC 60601 1	LGA	Inconceivable	
2.3 Electromagnetic hazards (according to IEC 60601-1-2)														
HA500	Electric fields hazards	Interruption of treatment (continuation not possible, further treatment necessary)	critical	Electromagnetic interference between system components and external components	occasional		Performance of the emission and immunity test for the device according to IEC 60601-1 clause 17 and 60601-1-2.	IEC 60601-1:2005 17 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	EMC-Report of IEC 60601-1-2	SLG	Inconceivable	
HA510	Electric fields hazards	Complaint	minor	Missing table(s) in user manual	hardly conceivable		Note in user manual: The user manual contains EMC table(s) for emissions and immunity according to DIN EN 60601-1-2	IEC 60601-1:2005 17 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	User Manual	ZI	Inconceivable	
HA520	Electricity (Application part)	Hazardous voltage in touchable	critical	Induction of energy interference of isolation by electric or magnetic fields	inconceivable		no measure necessary because risk is acceptable	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007				Inconceivable	
HA530	Electricity (Power supply)	Flash over to signal lines	critical	electromagnetic effects Harmed isolation	inconceivable		no measure necessary because risk is acceptable	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007				inconceivable	
HA540	Electricity (Source of supply)	Increasing ground leakage current	critical	Interference of isolation by electric/magnetic fields	inconceivable		no measure necessary because risk is acceptable	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007				inconceivable	
HA550	Emitted heat (Internal cooling system)	Overheating	critical	Interference of function by electric or magnetic fields	inconceivable		no measure necessary because risk is acceptable	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007				inconceivable	
HA560	Non-ionizing radiation (Internal electronics including electronic shutter)	Unintended change of adjusted parameters	critical	Induction of energy Interference of function by electric or magnetic fields	hardly conceivable		Shielding of the device, Filtering of the leads, ESD assembly	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007	Ja	EMC-Report of IEC 60601-1-2	SLG	Inconceivable	
HA570	Magnetic fields (Magnetic sensors)	Generally malfunctions	critical	Interference of function by external magnetic fields	hardly conceivable		System design according to IEC 60601-1-2	IEC 60601-1-2:2007	IEC 60601-1-2:2007 / EN 60601-1-2:2007	Ja	EMC-Report of IEC 60601-1-2	SLG	Inconceivable	

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Hazardanalyse / Risk Analysis						Risikobewertung / Risk Control											
Nr. Gefährdung No. Hazard	Med./techn. Auswirkg. Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Erstbewertung First evaluation	Maßnahme(n) Counter measure(s)	Angewandte Normen Applied Standards	Sicherheit Safety				Umsetzung der Maßnahmen Counter Measures Implementation		Restrisikobewertung Residual Risk Evaluation			
								Konstruktiv/Design	Schutzmaß./Safeguard	Information	Wirksamkeit geprüft? Effectiveness tested?	Umsetzung geprüft? Implementation checked?	Nachweisdokument(e) Record(s)	Verantwortlicher Responsible person	A	Restrisikobewertung Evaluation of residual risk	Risiko-Nutzen-Analyse Risk/Benefit Analysis
HA680	VIS radiation energy is too high	General Injury (moderate)	critical	Incorrect output	often	The safety limits (Fluence, pulse duration, repetition frequency, power) are supervised by hard- and software measures. - SW-watchdog - Checksum against file error - Ready-Standby procedure - Master-Slave communication (mutual monitoring, independent systems, plausibility of parameters) - Acoustical signal if switch to ready mode - Activation of laser is possible two seconds after activation of ready state - laser emission Indicator on GUI - Laser warning lamp - key switch - Door lock - Laser stop button - Redundant monitoring of output power	IEC 60601-1:2005 12.4.4 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja	Ja	Ja	Verification Report(s)	MMo	inconceivable		
HA690	No ionising radiation – laser radiation	Eye damage / injury, Injuries by fright	critical	Automatic release	occasional	- Stand by / ready procedure - Door interlock - Redundant release for software and hardware - Fiber detection	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		
HA700	No ionising radiation – laser radiation	Eye damage / injury, Injuries by fright, Injury of third persons	critical	Footswitch released unintentionally (laser radiation)	occasional	- Stand by / ready procedure - Door interlock - Redundant release for software and hardware - Fiber detection - footswitch with cover - Deactivation of laser in periods of inactivity (3min)	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		
HA710	No ionising radiation – laser radiation	Eye damage / injury, Injuries by fright, Injury of third persons	critical	No interruption of laser radiation in case of foot-switch malfunction and/or short circuit in signal lines)	occasional	- Redundant release for software and hardware (Two independent contacts) - Laser stop button	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		
HA720	No ionising radiation – laser radiation	Eye damage / injury, Injuries by fright, Injury of third persons	critical	Unintended release of laser radiation at ready activation be-cause of active footswitch	occasional	No release of laser radiation at ready activation if footswitch is activated.	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		
HA730	No ionising radiation – laser radiation	Eye damage / injury, Injury of third persons	critical	People in the room without safety goggles during laser or pulse light treatment	occasional	Acoustical and optical display of emission. Eyelid closing time at ready activation (1-Acoustical signal after ready state is activated. 2- Activation of laser is possible two seconds after activation of ready state) Note in user manual about using of goggles during laser treatment.	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis User manual	NJ ZI	inconceivable		
HA740	No ionising radiation – laser radiation	Eye damage / injury, Injury of third persons	critical	Fiber is not connected	occasional	1. No fiber is connected before treatment is started: --> The ready mode is not available and can not be activated 2. Disconnection of fiber during treatment: --> Error message is displayed. Switch to standby mode	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		
HA750	No ionising radiation – laser radiation	Eye damage / injury, Injury of third persons	critical	Access of third persons during laser treatment	occasional	- Door interlock - Laser warning lamp	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable		

Risikoanalyse / Risk Analysis										Risikobewertung / Risk Control																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																														
No.																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																																								

Risikoanalyse / Risk Analysis										Risikobewertung / Risk Control						
Nr. Gefährdung No. Hazard	Med./techn. Auswirkg. Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Erkennung First evaluation	Maßnahme(n) Counter measure(s)	Angewandte Normen Applied Standards	Sicherheits Safety		Umsetzung der Maßnahmen Counter Measures Implementation		Residualrisikobewertung Residual Risk Evaluation				
								Konstruktiv/Design Information	Schutzmaßn./Safeguard Information	Wirksamkeit geprüft? Effectiveness tested?	Umsetzung geprüft/ Implementation checked?	Nachweisdokument(e) Record(s)	Verantwortlicher/ Responsible person	A	Restrisikobewertung Evaluation of residual risk	Risiko-Nutzen-Analyse Risk/Benefit Analysis
HA880	Excessive temperatures in equipment	Burns (minor)	minor	[SFC] Failure of fan	occasional	Using of an approved fan: Documentation in list of critical components Measurement of the temperature in the device and turn off overheated components	IEC 60601-1:2005 4.7 [TRF] MDD 93/42/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja			LoCC FB 7331 LoCC QuadroPRO V01 MMo121101-01 SW-Verification Report	MMo NJ	Inconceivable		
HA890	Hazard by thermal energy	Burns (minor)	minor	[SFC] Omission/wrong dimensioning of fuses or over current releases	occasional	Device contains fuses 5x 20 T6,3A H 250V~ Documentation in list of critical components Labeling Note in User manual	IEC 60601-1:2005 8.11.5 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			LoCC FB 7331 LoCC QuadroPRO V01 MMo121101-01 Labeling Drawing 149099020 Labeling QuadroStarPRO CL 7357 Geräteschilder QuadroPro V01 IS120920-01 User Manual	MMo ZI	Inconceivable		
HA900	Excessive temperatures in equipment	Burns (minor)	minor	Temperature of ME-Equipment is too high	often	In normal use under worst case conditions (including the maximum ambient operating temperature) no parts do exceed temperatures given in IEC 60601-1 clause 11.1 table 23 and table 24. Allowable maximum temperatures which are touchable (table 23): - User Interface (1-10s; glass): Max. temp. 66°C - Handpiece (1min < t; metal): Max. temp. 48°C Allowable maximum temperatures which have skin contact (table 24) - Spacer (t > 10min; metal): Max. temp. 43°C Testing according to clause 11.1.3 for device	IEC 60601-1:2005 11.1.1 [TRF] IEC 60601-1:2005 11.1.3 [TRF] MDD 93/42/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Verification Report(s)	MMo	Inconceivable		
HA910	Excessive temperatures in equipment	Burns (minor)	critical	[SFC] Fire enclosures of ME Equipment doesn't withstand fire.	often	Insulated wire within the fire enclosure [CC] shall have a flammability classification equivalent FV-2, or better (according to the appropriate parts of the IEC 60695) [CC] - Connectors printed circuit boards and insulating material on which components are mounted [CC] shall have a flammability classification FV-2, or better, according to IEC 60695-11-10. [CC] - The fire enclosure meets the requirements of IEC60601-1, 11.3. - Critical components and fire enclosure parts are listed in "List of Critical Components - QuadroStarPro".	IEC 60601-1:2005 11.3 [TRF] MDD 93/42/EEC	IEC 60601-1:2005 / EN 60601-1:2006	Ja	Ja			MMo MMo	Inconceivable		
HA920						n.a. because there are no thermal cut-outs and over current releases with automatic resetting in the device	IEC 60601-1:2005 15.4.2.1a) [TRF]	IEC 60601-1:2005 / EN 60601-1:2006								
HA930						n.a. because there are no thermal cut-outs with a safety function that have to be reset by a soldering operation that can affect the operating value in the device.	IEC 60601-1:2005 15.4.2.1b) [TRF]	IEC 60601-1:2005 / EN 60601-1:2006								
HA940						n.a. because there are no thermostates in the device.	IEC 60601-1:2005 15.4.2.1c) [TRF]	IEC 60601-1:2005 / EN 60601-1:2006								

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Riskanalyse / Risk Analysis						Risikoherrschung / Risk Control											
Nr. Gefährdung No. Hazard	Med./techn. Auswirkg. Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Erstbewertung First evaluation	Maßnahme(n) Counter measure(s)	Angewandte Normen Applied Standards	Sicherheits- Safety				Umsetzung der Maßnahmen Counter Measures Implementation		Restrisikobewertung Residual Risk Evaluation			Risiko-Nutzen-Analyse Risk/Benefit Analysis
								Konstruktiv-Design Schutzmaßn./Safeguard Information	Wirksamkeit geprüft? Effectiveness tested?	Umsetzung geprüft? Implementation checked?	Nachweisdokument(e) Record(s)	Verantwortlicher Responsible person	A	Restrisikobewertung Evaluation of residual risk			
5.1.1 Hazards from falling caused by instability																	
HA1050						n.a. because QuadroStarPro is no mobile device.	IEC 60601-1:2005 9.4.2.4.3 [TRF] MDD 93/42/EEC	Headline IEC 60601-1:2005 / EN 60601-1:2006									
HA1060						n.a. because QuadroStarPro is no mobile device.	IEC 60601-1:2005 15.3.5 [TRF] MDD 93/42/EEC	IEC 60601-1:2005 / EN 60601-1:2006									
5.1.2 Suspended masses hazards																	
HA1070	Mechanical energy hazards	Injury as a result of excessive force (minor)	minor	Unintended movement of components.	occasional	Securing and fixing of components with latching connectors, cable ties and screws.	IEC 60601-1:2005 8.10.1 [TRF] MDD 93/42/EEC	Headline IEC 60601-1:2005 / EN 60601-1:2006	Ja			assembly drawings	KI	Inconceivable			
HA1080						n.a. because there is no support system.	IEC 60601-1:2005 9.8.1 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1090						n.a. because there is no support system.	IEC 60601-1:2005 9.8.2 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1100						n.a. because there is no support or suspension system.	IEC 60601-1:2005 9.8.3.1 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1110						n.a. because there is no support or suspension system.	IEC 60601-1:2005 9.8.3.2 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1120						n.a. because there is no support or suspension system.	IEC 60601-1:2005 9.8.4.1a) [TRF] IEC 60601-1:2005 9.8.4.1b) [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1130						n.a. because there is no support or suspension system.	IEC 60601-1:2005 9.8.4.3 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1140						n.a. because there is no support or suspension system.	IEC 60601-1:2005 9.8.5 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006									
HA1150	Mechanical energy hazards	Injury as a result of excessive force (minor)	minor	[SFC] Insufficient mechanical strength	often	QuadroStarPro have adequate mechanical strength and will not result in an unacceptable risk due to push test according to IEC 60601-1, item 15.3.2: - External parts of the enclosure shall withstand a force of 250N +/- 10N for 5s applied over a circular plane surface 30 mm in diameter.	IEC 60601-1:2005 15.3.2 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Test Report of IEC 60601-1	LGA	Inconceivable			
HA1160	Mechanical energy hazards	Injury as a result of excessive force (minor)	minor	[SFC] Insufficient mechanical strength	often	QuadroStarPro have adequate mechanical strength and will not result in an unacceptable risk due to impact test according to ISO 60601-1, item 15.3.3: - a solid steel ball (55 mm in diameter, mass of 500 g +/-25g) is permitted to fall freely from a 1,3 m height once onto each relevant part -> housing of device - To test vertical surfaces, the steel ball is suspended by a cord and allowed to swing like a pendulum in order to apply a horizontal impact, dropping through a vertical distance of 1,3 m once against each relevant part. -> housing of device	IEC 60601-1:2005 15.3.3 [TRF]	IEC 60601-1:2005 / EN 60601-1:2006	Ja			Test Report of IEC 60601-1	LGA	Inconceivable			

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Produkt-Risikomanagement/Product Risk Management

Risikoanalyse / Risk Analysis						Maßnahmenherleitung / Risk Control										Restriktionsbewertung / Residual Risk Evaluation		
Nr. Gefährdung No. Hazard	Med./techn. Auswirkung Medical/techn. Effect	S	Ursache der Fehlfunktion Cause of the malfunction	A	Erfassung First evaluation	Maßnahme(n) Counter measure(s)	Angewandte Normen Applied Standards	Sicherheit Safety	Umsetzung der Maßnahmen Counter Measures Implementation	Verantwortlicher Responsible person	A	Restriktionsbewertung Evaluation of residual risk	Risiko-Nutzen-Analyse Risk/Benefit Analysis					
HA1550	Operating outside required environmental conditions - overheating of the device - deviation of laser parameters - over treatment	minor	- Heat sink temperature too high	occasional		Handpiece: - Redundant temperature sensors on the warm side Laser: - Redundant temperature sensors on the warm side If the temperature too high the device switch in standby state. No hazard for patient, user and third parties.	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable			
HA1560	Failure of the device/function because of function/assembly faults	minor	- failures of current sensors	occasional		- Monitoring of admissible values - Redundance (current driver and photodiode)	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable			
HA1570	Hazards by investigation in the FDA data-base	minor	- energy output too high	occasional		- Monitoring of admissible values - Redundance (current driver and photodiode)	EN 62304: 2006	Ja	Ja	Ja	Ja	Ja	SW-Verification Report(s) Traceability analysis	NJ	inconceivable			
HA1580	Operational hazards (Uncontrolled function)	minor	- Under- or overvoltage	occasional		Wide range power supply with input voltage range of 100-240V (+/- 10%)	MDD 93/42/EEC	Ja	Ja	Ja	Ja	Ja	Datasheet	MMo	inconceivable			
HA1590	Operational hazards	minor	- Damaging of parts and accessories due to multiple cleaning and disinfection.	occasional		Note in user manual about multiple cleaning and disinfection procedure for QuadroStarPro	IEC 60601-1:2005 11.6.6 [TRF] MDD 93/42/EEC	Ja	Ja	Ja	Ja	Ja	User Manual	ZI	inconceivable			
HA1600						n.a. because no sterilization necessary	IEC 60601-1:2005 11.6.7 [TRF] MDD 93/42/EEC	Ja	Ja	Ja	Ja	Ja						
HA1610	Operational hazards	minor	- Missing or incomprehensible formulated accompanying documents	occasional		Accompanying documents for QuadroStarPro consist of: - User Manual (UM) - Technical description as part of the User Manual The documentation is provided as hardcopy	IEC 60601-1:2005 7.9.1 [TRF] MDD 93/42/EEC	Ja	Ja	Ja	Ja	Ja	User Manual	ZI	inconceivable			
HA1620	Storage and operating out-side the specified ambient conditions	minor	- Insufficient packaging	occasional		Testing of packaging according to specified transport, storage and operating conditions	MDD 93/42/EEC	Ja	Ja	Ja	Ja	Ja	Verification Report Packaging	KI	inconceivable			
HA1630	Well-founded predictable misuse (non-intended use)	critical	- Human error	occasional		Intuitive design of QuadroStarPro. Testing of usability.	EN 60601-1-8: 2004	Ja	Ja	Ja	Ja	Ja	Usability-File		inconceivable			
HA1640	Insufficient warning of side effects (Wrong treatment parameters)	critical	- Not or badly skilled user	occasional		Description of user profile for intended user Training of personnel	EN 60601-1-8: 2004	Ja	Ja	Ja	Ja	Ja	Usability-File User manual Training protocol		inconceivable			
HA1650	Insufficient warning of contraindication (Misuse of the device)	minor	- Uninsufficient description of intended use	occasional		Description of intended use in user manual	EN 60601-1-8: 2004	Ja	Ja	Ja	Ja	Ja	User manual		inconceivable			
HA1660	Transport damage (Too high weight or bad weight spreading; Rough edges)	minor	- Too big, non-ergonomic designed and weighty device, packaging insufficient	occasional		Ergonomic and small design of the device. Lightweight (< 15kg) Package testing	EN 60601-1-8: 2004	Ja	Ja	Ja	Ja	Ja	Zeichnung User manual verification Report for packaging		inconceivable			
8 Labelling																		
8.1 Hazards related to installation instructions																		
HA1670						No hazard during primature unpacking of the QuadroStarPro.	IEC 60601-1:2005 7.2.17 [TRF]	Ja	Ja	Ja	Ja	Ja						
8.2 Hazards related to device labelling																		

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Form 7301

Risk Management Concept / Report acc. to DIN EN ISO 14971/ ISO 14971

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Index: C

Identification number: IS/12/07/24-02

Member code/year/month/day-sequential number

 Page 1 of 5
 FB 7301 Risk Management Concept
 QPro IS120724-02.doc

For:

Medical product: QuadroStar Pro

Item number: 1490 - 1496

Assessed construction level: --

Revision users manual: --

☐ Before R&D☒ Risk management concept☒ During R&D☐ Risk management report☐ After R&D (before serial production)☐ After production

Drawn up:

Vs.	Description of change	Date	Name	Sign
01	First edition	July 24, 2012	I. Streit	
02				

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Not applicable for risk management concept.

2. Responsibilities

Responsibility	Name	Qualification	Responsible in the project
overall responsibility acc. to DIN EN 60 ISO 14971	M. Möller	Dipl.-Ing. Biomedizintechnik (FH)	Safety management total device
electrical safety acc. to DIN EN 60 601-1, DIN EN 60 601-1-2 DIN EN 60 601-2-22	C. Büßer	Dipl.-Ing. Elektrotechnik (FH)	Electric safety EMC
software acc. to DIN EN 60 601-1, EN 62304	N. John	Dipl.- Ing. Technische Informatik (FH)	Software
Usability acc. to EN 60601-1-6, EN 62366:2008	G. Zimmermann	Physicist	Usability, Application
laser safety acc. to DIN EN 825-1, DIN EN 60 601-2-22	I. Streit	Physicist, Dr. rer. nat.	Laser safety

All risk management documents are checked and released by the quality representative QMB.

3. Criteria for defensibleness of risks

The following matrix has to be applied for risk assessment of all detected hazards.

Occurrence probability

always				
often				
occasional				
hardly conceivable				
inconceivable				
	dispensable	minor	critical	disastrous

Extent of damage

Accepting criteria


	not passable	ALARP	passable
--	--------------	-------	----------

ALARP = As Low As Reasonable Practicable

Damages by contact with supply voltage and eye damages caused by laser radiation are eminently critical at devices distributed by Asclepion Laser Technologies. These risks have to be looked particularly carefully therefore.

Form 7301	Risk Management Concept / Report acc. to DIN EN ISO 14971/ ISO 14971 (In Jewells gültiger Ausgabe entsprechend CL7321 Relevante Sicherheitsnormen)	 Asclepion Laser Technologies
Index: C	Identification number: IS/12/07/24-02 Member code/year/month/day-sequential number	Page 4 of 5 FB 7301 Risk Management Concept QPro IS120724-02.doc

Extent of damage and Occurrence probability are selected in the following way:

Form 7301	Risk Management Concept / Report acc. to DIN EN ISO 14971/ ISO 14971 (in jeweils gültiger Ausgabe entsprechend CL7321 Relevante Sicherheitsnormen)	
Index: C	Identification number: IS/12/07/24-02 Member code/year/month/day-sequential number	Page 5 of 5 FB 7301 Risk Management Concept QPro IS120724-02.doc

4. Risk management process

The process of the risk management during the product development phase is carried out in accordance with the process instruction PA 731 "Design and development planning" and the check list CL 7311 "Development plan", or else in accordance with process instruction PA 737 "Design changes".

5. Intended use

Not applicable for risk management concept.

6. Risk analysis and risk assessment

The risk analysis generally consists of the components:

- general risk analysis according to the Medical Device Directive **MDD 93/42/EWG:1993-06+ M1+M2+M3+M4+M5**
- risk analysis according to **EN 60601-1:2006** for general safety
- risk analysis according to **EN 60601-1:2006** and **EN 62304:2006** for PEMS (safety of software).
- risk analysis according to **EN 60 601-1-2:2007** for the determination of the essential performance features for the specification of the parameters in the interference immunity tests
- proof of usability according to **EN 60601-1-6:2007+CENELEC-Cor.:2010** and **EN 62366:2008** for the use of medical device and contributed factors

7. Description of the applied risk decreasing measures

Not applicable for risk management concept.

8. Risk control and risk verification

The verification and validation of the risk measures was carried out according to the process instructions PA 735 "Design Verification" and PA 736. "Design validation".

Form 7301

Risk Management Concept / Report acc. to

EN ISO 14971: 2009

Asclepion
Laser Technologies

Index: E

Identification number:

MMo/13/10/01-02

Member code/year/month/day-sequential number

Page 1 of 29

FB 7301 RMR QuPRO
V02_MM0131001-02.doc

For:

Medical product:

QuadroStarPRO

☐ Before R&D

Item number:

1490 - 1495

☐ Risk management concept☐ During R&D

Assessed construction level:

Pre series

☒ Risk management report☒ After R&D (before serial production)

Revision users manual


1490 99 001_Vs2-10.Febr14-GA_QuPRO-DE

☐ After production

Drawn up:

Vs.	Description of change	Date	Name	Sign
01	First edition	12.11.2012	M. Möller	Sgn. Möller
02	No pressure sensor, Bypass for pressure limitation (8.1.1.15.1; 8.1.8) Risk management acc. EN ISO 14971:2009	20.02.2014	M. Möller	

Checked and released:

	Division	Date	Name	Sign
Checked by:	R&D	20.02.2014	N. John	
Released by:	R&D	Feb. 21, 2014	I. Streit	

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Datum

Name

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Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 2 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

0. Contents

Description of change	1	8.1.1.12.2 TEC system	17
0. Contents	2	8.1.1.13 Monitoring of hard and software status	17
1. Summary / Zusammenfassung	3	8.1.1.14 Efficiency of laser power output (current/power adjustment)	18
2. Description of the device / Beschreibung des Gerätes	6	8.1.1.15 Electronic shutter	18
2.1 Brief description / Kurzbeschreibung	6	8.1.1.16 Monitoring of current and photodiode value	19
3. Responsibilities	8	8.1.1.17 Fiber recognize	21
4. Criteria for defensibility of risks	8	8.1.1.18 Fiber detection	21
5. Risk management process	8	8.1.1.19 Safety Systems of the Scanner unit	21
6. Intended use	8	8.1.2 Protection against Electrical hazards	22
7. Risk analysis and risk assessment	9	8.1.3 Protection against electromagnetic phenomena	22
8. Description of the applied risk decreasing measures	9	8.1.4 Protection against laser radiation	23
8.1 Risks	9	8.1.5 Protection against unwanted laser emission	23
8.1.1 Protection against Software failures	10	8.1.6 Protection against wrong use	23
8.1.1.1 Watch-dog system	10	8.1.7 Protection against hazardous or wrong output	23
8.1.1.2 Monitoring clock frequency	12	8.1.8 Failure of mechanical, resistance to heat, expelled parts	26
8.1.1.3 Checksums over files	12	8.1.9 moving parts	27
8.1.1.4 Master-Slave communication	13	8.1.10 Movement	27
8.1.1.5 Laser emergency stop button, key switch and power switch	13	8.1.11 Safety signs	27
8.1.1.6 Ready-Standby circuitry	15	8.1.12 Cleaning and disinfection, Sterilization and Compatibility	28
8.1.1.7 Foot button with cover and two contact pairs	15	8.1.13 Accuracy and Arrangement of controls and indicators	28
8.1.1.8 Door interlock	15	8.1.14 Acoustic energy	28
8.1.1.9 External laser warning lamp	16	8.1.15 Batteries	28
8.1.1.10 Optical and acoustic laser emission indication	16	9. Risk from inadequate controlled processes in design, manufacture and maintenance	29
8.1.1.11 Eyelid closing time	16	10. Risk control and risk verification	29
8.1.1.12 Temperature measurement (device)	16		
8.1.1.12.1 Skin Cooling system (Water flow, Temperature monitoring)			

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02	Page 3 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

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Remarks: The risk management of the device was carried out according to EN ISO 14971:2009. At present, an analysis of the implementation of EN ISO 14971:2012 is performed. For this purpose, it is planned to adapt the PA 730 to perform a GAP analysis to the previous made risk management and adapt the risk analysis.

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For the "QuadroStarPRO" are several essential performance identified. These performances should be checked by inspection, and if necessary, by functional test. The essential performances are:

- Laser emission wavelength / Laserstrahlung Wellenlänge:
 - Laser emission wave length: 532nm (green) or 577nm (yellow) or 940nm or 980nm (infrared)
 - Laser energy and pulse length/ Laserenergie und Pulslänge
 - Laser emission power: 8W @ 532nm (green) or 5W @ 577nm (yellow) or 25W @ 940nm or 980nm (infrared)
 - Laser pulse length: 1ms up to 95sec, CW
 - Repetition rate: 0,5 – 20Hz
 - Laser class: 4
- *for fiber with core diameters smaller than 600µm, the maximum power is limited automatically

The output parameters and the safety systems of the device are monitored and controlled by a microcontroller system. Equivalent voltage information will be evaluated and compared with desired value.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 4 of 29 FD 0001 FORM 03-07-10 VOL. 0001 FORM 03-07-10
		Member code/year/month/day-sequential number	

See	Limit	Result
Protection against hazardous or wrong output	Nominal energy < $\pm 20\%$	warning sign will be shown at the display
	Nominal energy < $\pm 50\%$	pulse generation will be stopped immediately and an error information will be shown
	Nominal pulse duration < $\pm 20\%$	switch off the laser
	failure tolerance time of the system > watchdog max. 200 ms	If the laser emission will interrupted after this time no unacceptable injury is possible.

- Isolation Applied Part Type B / Isolation Anwendungsteil Typ B
The applied part is a type B part. The appliance is designed for use in patient environment. (look at insulation diagram)

See	Result
Pkt. 8.1.2 Protection against Electrical hazards	- an reinforced insulation of input line voltage (life part) from the safety low voltage supply inside the device and the applied part with 2 MOPP - Basic Insulation (1 MOPP) between safety low voltage and grounded applied part - On the backside metallic enclosure with potential earth, basic insulation (1MOPP) life part to enclosure - Grounded applied part

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Laser class: 4
 Power supply: 100 - 240VAC, 50/60 Hz
 Fluids: It is not intended to use fluids together with the device
 Ambient conditions: Device is mainly used in doctor's office, but in OR as well
 Handpiece and scanner are used in doctor's office – treatments of intact skin (distance holder of the handpiece can be disinfected or sterilized)
 Fibers can be used in OR for coagulation and vaporization of tissue under sterile conditions also – fibers are sterile
 User of device: practising physician

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 5 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

1.3 Assessment of the residual risk / Bewertung des Restrisikos

Verbal description

Using the product "QuadroStarPRO" for its intended use acc. to the user Manual does not include any critical risks to patients, user and innocent third parties.

By the laid out specifications, built-in security measures, Quality management actions and quality-oriented production activities the remaining risks due to technical factors of the product are acceptable.

Through maintenance and safety checks, the performance, specification and safety of the product can be guaranteed for the intended service life.

In the risk analysis and assessment of clinical the residual risks are assessed to be acceptable, all measures set out in the risk management process were verified positive.

The process defined in the QMH (PA 730 Risk Management) to control the risk is appropriate and sufficient.

The documented risk management process provides acceptance and post market monitoring of all risks for the product "QuadroStarPRO".

All at this moment known hazards were assessed.

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☒ yes / ja ☐ no / nein ☐ n/a

All measures were implemented and verified.

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☒ yes / ja ☐ no / nein ☐ n/a

The medical benefit is greater than the residual risk.

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Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 6 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

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
The QuadroStarPRO is a medical therapeutic laser of Class 4. Different wavelengths are available:

- 532nm (green) or
- 577nm (yellow) or
- 940nm or 980nm (infrared)

The device is equipped with a power switch for 100-240VAC 50/60 Hz for switching it on and off. For its intended use the device must be activated by a key switch. A laser emergency button on the front cover is provided for the emergency case. The device is controlled by a master-slave software system. The master process is a graphical user interface running at a Marvell PXA-300 CPU powered by Windows CE operating system. The Slave is a separate process running at a Microchip PIC18F46K22 Microcontroller. The communication between master and slave is done via secured inter-process communication. The slave control unit is for the electronics of the laser. For a comfortable human-machine interaction a 5.7" color display with 9 capacitive touch buttons around is used at the device. The Master software is created by means by Microsoft Visual Studio 2008 and the Slave Software is created by MPLAB X IDE.

The Master CPU is plugged onto a DBP1 (Display Backplane) board developed by Asclepion Laser Technologies GmbH (see corresponding part specification).

The complete electronic system works with low voltage. A foot switch is used to release the laser emission.

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Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 8 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

3. Responsibilities

See chapter 2 in FB 7301 Risk Management Concept QPro IS120724-02

4. Criteria for defensibility of risks

See chapter 3 in FB 7301 Risk Management Concept QPro IS120724-02

5. Risk management process

See chapter 4 in FB 7301 Risk Management Concept QPro IS120724-02

6. Intended use

The QuadroStarPRO is a laser (Class 4 laser product) that has relatively long impulses up to a continuing emission at all varieties of wavelengths. Therefore the device can cause a selective thermal impact or damage to the tissue and hence can be used for various indications. The versions in the visible (green and yellow) and in the infrared range are nearly identical in their effect. For this reason the intended use and indications are given together for the specific wavelength. The intended use is described in more detail in chapter 1.3 and 6 of the User's Manual QuPRO (item 1490 99 002).

This device may only be applied by Users who are adequately qualified to do so, based on their professional training background, expert knowledge and practical experience. Treatment may not be performed by anyone other than a practicing physician or a well trained person like a nurse under the responsibility of a practicing physician. This device contains in particular a high power laser system. Laser knowledge becomes presupposed or must become acquired before the use of this device. Appropriate training courses to provide laser expert knowledge are also offered by Asclepion Laser Technologies GmbH. Only those users are permitted to use this device which have received on-the-job instructions and signed off for them in the Medical Device Log-book.

The device can be used in any private practice or clinic and easily relocated to the desired place of treatment. The QuadroStarPRO is designed to the generally accepted rules of technology and meets all relevant safety standards. The software is guided and can only be changed by authorized service personnel.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 9 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

7. Risk analysis and risk assessment

The risk analysis generally consists of the components:

- general risk analysis according to the **MDD 93/42/EWG:1993-06+ M1+M2+M3+M4+M5**
- risk analysis according to **EN 60601-1:2006** for General requirements for safety
- risk analysis according to **EN 62304:2006** for PEMS (safety of software)
- risk analysis according to **EN 60601-1-2:2007** for the determination of the essential performance features for the specification of the parameters in the interference immunity tests
- risk analysis according to **EN 60601-1-6:2007+CENELEC-Cor.:2010, EN 62366:2008** for the determination of the usability
- risk analysis according to **EN 60601-2-22:1996, EN 60825-1:2007** for Particular requirements for the safety of therapeutic laser

The QuadroStarPRO was checked with a type test. Other test or procedures have not been established.

8. Description of the applied risk decreasing measures

8.1 Risks

Like for all devices, risks are to be taken into consideration when using this device.

The risks for the patient, the user, the service personnel and bystanders include

- Failure of security managements (by Software)
- electrical shocks
- risks by electromagnetic phenomena
- damage to the eyes by uncontrolled emission or faulty output values of laser light
- failure of mechanical integrity

The QuadroStarPRO has been developed, constructed and manufactured in such a way that it does not cause any foreseeable risks for the clients or for persons and objects being in its environment, provided the module is used as intended according to the User Manual.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 10 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

8.1.1 Protection against Software failures

8.1.1.1 Watch-dog system

A computer system can only partly monitor itself. A watchdog system cares about the required reliability. It is an external unit operating independently of the computer. Its task is to switch the PEMS into a reliable mode in case of missing control signals. Control signals will be missed, if:

- the software is in an endless loop
- any content of the memory is damaged
- the peripheral components of the CPU have a technical defect
- the CPU doesn't process the prescribed program path due to external or internal influences
- CPU and/or system errors cause the stop of the CPU
- A single and/or systematic software error has been detected.

The software starts automatically after power on. The time of the watch dog unit is checked several times by the software. If the time is in valid range the software starts triggering the watch dog. As long as the unit is triggered correctly, the output enable signal is enabled (→valid state). If the software is not able to trigger the watch dog during program execution, all outputs will be switched to safe state due to output enable signal is disabled (→invalid state). Once the refreshing of the watch dog failed it can't be restarted during program execution. The state of output enable signal is always read back by the Controller and checked before a new watch dog trigger is generated. If the signal is disabled (invalid) no watchdog triggers will be generated. The device must be restarted by power-off/on.

Watchdog 1

The first watchdog is an internal PIC watchdog. All outputs will be cleared if no signal occurs more than 4s after µC start.

Watchdog 2

The watchdog is cleared when any of the following events occur: a SLEEP or CLRWDT instruction is executed.

The external watchdog is trigger by "KeepAlive-Flag". This event flag is toggled by software main loop and any other functions (for example Communication). The longest time to refresh watchdog is 100ms trigger through main loop. If Main loop stopped and timer increase more than 400ms watchdog set device OFF.

Form 7301	Risk Management Concept / Report acc. to		
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 11 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	


Figure 1

Malfunction 1: The interrupt fails

In this case, the rectangular signal does not exist any longer. The watchdog trigger will not be started again. The watchdog timer goes down.

- All output signals will be reset.
- The control system is safe and stays on.
- The device turns off by itself after 30 seconds.

Malfunction 2: The software crashes and the interrupt remain active.

Due to the missing refresh flag the refresh time is not reset to 110ms in this case and runs down to zero. Therefore, the gate for the rectangular signal will be closed and another will not be generated. Consequently, there is no watchdog refresh and consequently the watchdog timer goes down.

- All output signals will be reset.
- The control system is safe and stays on.
- The device shows an error message and turns off by itself after 30 seconds.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 12 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

Maximal time:

Watchdog refresh main all 100ms in loop. Whenever time is more than 400ms device shutdown.

Minimal time:

→ Whenever the device is switched on, the watchdog timer will be checked several times if the watchdog timer in range: $t \leq 400\text{ms}$

8.1.1.2 Monitoring clock frequency

Controller timers are used for system time base and triggering of several events (for e.g. watch dog interrupt). Accuracy of timers depends on stable Controller frequency. A Controller independent external is UART-signal to communicate to master. If clock frequency deviation exceeds one percent ready mode is blocked by an error message.

8.1.1.3 Checksums over files

Procedure: In this procedure, the contents of the memory/the data are reduced by the arithmetic addition of all bytes to one sum with a number of digits (32bit) depending on the memory capacity. Carries are also taken into consideration. In this sophisticated and quick test procedure it must be ensured that faulty information changes in the memory take place only towards one direction, i.e. always from "0" to "1" or always from "1" to "0" only! This feature cannot be guaranteed for the system. That means that bit errors can cancel each other out and are not detected. The probability of bit errors cancelling each other out is seen as very low and the procedure is considered appropriate and sufficiently reliable.

The CPU module is provided with a flash memory for data archiving purposes. This memory is configured as *FlashDisk* and can be addressed via a flash file system. Thus, data can be accessed via files as on usual hard disk drives by using the functionality of the operating system. To identify read errors and/or memory errors all files created by Asclepion Laser Technologies GmbH are provided with a check sum.

At the end of the development, the software developer stores all files with their corresponding check sums (stored at the end of each file). Whenever a file is accessed, the software calculates the belonging check sum and compares it to the check sum stored at the end of the file. If both sums are identical the data is valid and can be read and used.

...the software calculates the belonging check sum and compares it to the check sum stored at the end of the file. If both sums are identical the data is valid and can be read and used.

...the developer stores all files with their corresponding check sums (stored at the end of each file). Whenever a file is stored at the end of the file. If both sums are identical the

Form 7301	Risk Management Concept / Report acc. to		 Asclepion Laser Technologies
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 13 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

The user software is only started if check sums of all needed files are valid. The user software (taurus.exe) tests the integrity of used configuration files (INI files) during runtime whenever they are required. In case of any error default data is used and an error message is shown. The ready mode activation of the belonging module is blocked.

8.1.1.4 Master-Slave communication

User-device interaction is realized using buttons next to display. This interaction is controlled by the master software which provides the graphical user interface. The master is running on separate CPU. The slave runs on μ -Controller. The master can't access the device hardware directly - this is only possible in the slave process. Therefore master and slave are connected using a secured inter-process communication.

If the slave process crashes the device is secured by the watch dog. If the master process crashes the slave process switches the device to safe mode and tries to re-establish the communication to master. If this attempt fails the slave turns off the device. The communication is bidirectional. Every transmission is secured by a CRC for the transmission data. The master checks the response of the slave after any command for being correct. The slave expects a request from the master within a defined time period. For any erroneous transmission the transmission cycle is repeated. If erroneous transmission occurs continuously on repeat the slave turns off the device.

Any malfunction detected by slave is signalized to master (state bit). Thereafter the master requests and displays the error information.

If the ready button at the display is pressed by user. The master calculates the necessary treatment control values and monitors the emission of energy. The master transmits data with all treatment parameters to slave. Slave configures hardware and sends set parameter back to master. Master check set data read back from slave. If different parameters determined master generate a message and stop ready procedure.

8.1.1.5 Laser emergency stop button, key switch and power switch

The laser emergency button and key switch are in line of the power supply. Both turn on the voltage from the power supply relay K1 at the SOC PCB. If the emergency stop button is pressed, the power supply is OFF, the system is safe.

The key switch is an on-off-switch that is evaluated by the PEMS. If the switch is in the position OFF, the PEMS switches itself. It prevents the device from being operated by unauthorized persons.

The power switch disconnects the device from the line voltage (separates all connections). It prevents the STANDBY mode, i.e. the device can only be activated by the key switch if the power switch is in the ON position.

Form 7301		Risk Management Concept / Report acc. to EN ISO 14971: 2009	Asclepion Laser Technologies
Index: E		Identification number: MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 14 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc


Figure 2

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 15 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

8.1.1.6 Ready-Standby circuitry

After power-on the device will be initialized by hard- and software. The device is always in standby mode after software start. After changing a module it is always in standby mode, too.

The user must explicitly switch on READY mode for the selected module by using the button at the display. The READY mode is clearly indicated at the display. All treatment parameters are positioned in the field of vision near to the READY button. The user is responsible for assessing and checking them visually when activating the READY mode.

Treatment control signals of the selected module are only enabled if READY mode was activated successfully. If STANDBY mode was activated, all treatment control signals are disabled. After a period of inactivity during ready mode (foot switch not pressed for 3 minutes) the selected module is switched to standby mode automatically.

In case of any error the device will be switched to standby mode and error information is shown on display.

8.1.1.7 Foot button with cover and two contact pairs

In order to avoid unintentional activation of the foot switch by the user, third persons or falling down objects, the foot switch has been equipped with a protective clamp.

Due to redundancy, the foot switch is equipped with two separated contact pairs. Both signals are evaluated separately by the PEMS. The system will only recognize the footswitch as being pressed when both contact pairs are activated simultaneously.

Whenever ready mode is activated, the foot switch will be checked for one of both contacts being active. In this case an error is generated and device remains in standby mode.

8.1.1.8 Door interlock

The device is provided with an optional door interlock switch. This contact designed as a normally open contact is to be connected to the entrance(s). Activation of READY mode is only possible, if this contact is closed. If the contact is opens during READY mode, the laser emission is stopped immediately, STANDBY mode is activated and a warning message is displayed. Therefore, third persons will be protected when entering the laser room without wearing laser safety goggles.

- This check is performed when switching into the ready mode and continuously during READY mode.
- If the interlock contact is open, the system switches to the STANDBY mode and outputs an error message.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 16 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

8.1.1.9 External laser warning lamp

In order to prevent persons from entering the laser treatment room, the device offers a potential-free contact for signaling the READY mode. This signal can be used in front of the entrance to provide an optical and/or acoustic indication of potential danger in the laser treatment room.

8.1.1.10 Optical and acoustic laser emission indication

Whenever the device changes to READY mode, an acoustic signal can be heard. During the laser emission, an acoustic/optical signal is output synchronously with the pulse.

8.1.1.11 Eyelid closing time

According to valid regulations, all persons in the laser room must use the required laser-protective eye-wear. If one person does not wear a protective eye-wear, the change to the READY mode will be signaled optically and acoustically. During this time, the foot switch is not evaluated for two seconds (eyelid-closing time).

8.1.1.12 Temperature measurement (device)

The software uses multi-tasking. The temperature management is divided into three areas:

- Device temperature
- TEC temperature
- Skincooler temperature

The fans for cooling the device work permanently. The fans for cooling the heatsink are controlled permanently and work depending on the heatsink temperature. The status request from master calls device temperature and set device shut down if it's over.

The Temperature of the TEC modules is controlled by TEC module and Skin cooler temperature is controlled by Skin cooler module separate. The separate temperature modules (TEC, Skin) report the temperature state to the slave system. The slave evaluates the temperature state and set the device in a safe condition (e.g. standby), if the state is out of range. The status request from master calls temperature module state. If the state is out of range, a failure notice displayed.

Form 7301	Risk Management Concept / Report acc. to		 Asclepion Laser Technologies
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 17 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

8.1.1.12.1 Skin Cooling system (Water flow, Temperature monitoring)

Skin Cooling system works in a separate module. It communicates with the slave system. The cooling system itself works independently from slave system and reports the actual state to the slave system. The cooling system module can be controlled by the slave system.

Cooler continuously checks itself for:

- main water temperature is in valid range
- handpiece cold side temperature is in valid range
- state for coolant flow (check of temperature)

In case of any error cooler is disabled and an error message is shown.

8.1.1.12.2 TEC system

TEC Cooling system works in a separate module. It controls the Laser temperature and is configured redundantly. The TEC cooling system communicates with slave system. The cooling system itself works independently from slave system and reports the actual state to the slave system. The cooling system module can be controlled by the slave system.

Controller continuously checks:

- cold side temperature
 - controlled in ready or turn to ready
 - out of range → set device in standby and show error
- warm side temperature
 - controlled in ready & standby
 - out of range → set TEC and pump OFF, show error

8.1.1.13 Monitoring of hard and software status

Master generate loop to call slave state every 150ms. Master interface read received slave data and parse informations. The active laser module continuously performs all necessary checks for digital signals being valid. In cases of any malfunction information will be displayed. For worst case the watch dog is not refreshed and locks all outputs and set device OFF.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 18 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

8.1.1.14 Efficiency of laser power output (current/power adjustment)

The hardware to provide the laser diode current exists at the device inside. The components of this hardware are different from device to device. Therefore this hardware need an adjustment. For this purpose a current and power adjustment is necessary. The current/power adjustment will be stored (non-volatile) in the EEPROM of the slave controller. At this memory the current and power adjustment are stored in a set of points (drive value of the diode driver to laser diode current and to the power laser output). These points are predefined by a special service software for authorized service technician only.

The adjustment procedure adjusts a defined drive in CW Mode. The drive value is predefined by the software. The service technician measures the output power with a calibrated power meter and enters the measured power data into the service software. With entering the measured power the laser diode current will be automatically stored. If the procedure is complete the adjustment points will be calculated and stored in the controller.

During runtime in user mode this adjustment points will be used to calculate the drive in dependence of the adjusted output power.

8.1.1.15 Electronic shutter

To generate the laser radiation, the diode driver unit must be active, a control voltage (drive - PWM) must be applied and two trigger signals (enable, inhibit) must be produced. These control signals mainly depend on the slave software.

To increase the safety of the system, the laser has been provided with a separate electronic shutter system that operates independently of diode driver module. It is an electronic circuit (PSU PCB) that disconnects the 12V power supply of the diode driver.

Operating regime:

The right startup sequence consists of

1. Set Enable (ACTIVE HIGH): enable power supply and the power stage of the driver, bias current will flow
2. Set PWM level (PWM): set drive level
3. Pulsing the inhibit line (ACTIVE LOW): enables the PWM level, blocks current above the low bias current

State	PowerSupply (PSU PCB)	Enable	Inhibit	Current
KEYSWITCH OFF	OFF	OFF	OFF	ZERO
START	OFF	OFF	OFF	ZERO
STANDBY	OFF	OFF	OFF	ZERO

STANDBY	OFF	OFF	OFF	ZERO
READY	ON	OFF	ON	ZERO
OPERATE (Footswitch activated)	ON	ON	OFF	BIAS
		ON	ON/PULSE	Whatever set by PWM

Form 7301	Risk Management Concept / Report acc. to		
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 19 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

8.1.1.16 Monitoring of current and photodiode value

Slave checks

The drive value is set by Slave Controller and is proportional to:

- the current of the diode driver unit
- the output power of the laser module
- the photodiode value of the laser module

Diode Current (Slave)

- Every emitted pulse current, is measured asynchronously to the pulse interrupt in any of the following states (READY, OPERATE, SCAN) where the current driver is enabled
- the sampled value is compared with the programmed value
- the resulting error is compared with the range (12.5% of the programmed value)
- in case the error exceeds the range value the error count is increased. Otherwise it is decreased.
- in case the max. error value is exceeded the device is switched to standby and an error message is displayed

Photodiode (Slave)

- Every emitted pulse is measured asynchronously to the pulse interrupt in any of the following states (READY, OPERATE, SCAN) where the current driver is enabled
- Photodiode sample routine
 - o New value sampled every 100µs
 - o The sampled value passes through the software filter
 - o The value is compared with the programmed value
 - o The resulting error is compared with the range
 - o In case the error exceeds the range value the error count is increased. Otherwise it is decreased
 - o The error is compared with the max_error value (1000) to generate a variable error (error time 100ms)

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 20 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

Master checks

Photodiode (Master)

- Every received slave state in ready master compare value (set and get).
- Photodiode sample routine (READY)
 - o Read actual value
 - o Compare get and set value
 - o if the limit photodiode current 50% too large (positive or negative) the pulse generation will be stopped immediately and an error information will be shown
 - o if the limit photodiode current 20% too large (positive or negative) a warning sign will be shown at the display
 - if more than 20% set icon for over energy
 - if less than 20% set icon for under energy
- Photodiode sample routine (STANDBY)
 - o Read actual value
 - o If actual slave state STANDBY and photodiode current is more than 10 digits → show message and set device OFF

Diode current (Master)

- Every received slave state in ready master compare value (set and get).
- Photodiode sample routine (READY)
 - o Read actual value
 - o Compare get and set value
 - o if the limit photodiode current 50% too large (positive or negative) the pulse generation will be stopped immediately and an error information will be shown
 - o if the limit photodiode current 20% too large (positive or negative) a warning sign will be shown at the display
 - if more than 20% set icon for over energy
 - if less than 20% set icon for under energy

Form 7301	Risk Management Concept / Report acc. to		 Asclepion Laser Technologies
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 21 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

8.1.1.17 Fiber recognize

To ensure that the correct fiber is used, it is coded with 4-bits. This bit coding provides the identification of 15 hand pieces (0 is invalid). The fibers are assigned to the different modules with different power outputs. Without valid fiber no laser emission is possible. The energy density or power is always calculated, if a valid fiber is identified. The fiber recognition connection is checked permanently. Allows the fiber type different spot size, it can be selected in user menu.

8.1.1.18 Fiber detection

To ensure that the laser beam exits at the expected discharge opening of the applicator (handpiece), the fiber applicator must be screwed into the SMA coupling. To monitor the correct connection a micro switch is located on the fiber coupling system. The fiber detection connection is checked permanently. Without valid fiber connection no laser emission is possible.

8.1.1.19 Safety Systems of the Scanner unit

The scanner was designed to be compliant with requirements for medical systems. Two independent systems for communication an error situation to the host system can be used:

- 1) An error output signals any fault condition that requires the host system to shut down.
- 2) A ready output signal informs the host that the mirror is in position. This signal will not be output if an error condition exists.

The implemented safety systems are:

- The scanner module has closed loop mirror positioning, sensors continually keep track of the mirror positions and signals an error if the position is not as expected.
- Any major fault in the drivers, feedback loops or mechanical system will stop operation.
- The scanner module has an AutoPosition mode that uses the current mirror position to verify that the intended position has been reached. Only if the position is OK the scanner will signal that it is ready instead of relying on a fixed positioning time.
- Optical layout ensures that even if there is a mirror fault the optical beam will not pass through the scanner.
- A sensor on the mirror mount will stop laser operation if the temperature is too high. This can limit destructive damage in high power applications if the mirror fails or gets contaminated.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 22 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

- An error output from the scanner to the laser system signals if any fault condition is present. This error signal should stop laser output from the host system.
- Internal sensors stop the scanner and signals an error if the temperature is out of the safe operating range.

8.1.2 Protection against Electrical hazards

The device is designed as electrical class I equipment. The appliance is not a permanent installed device. Concerning leakage current the appliance fulfils the demands of EN 60601-1:2006. The maximum nominal line input is 240 Vac. Inside the device there is after the ac/dc power supplies only safety low voltage < 30 Vdc, reinforced insulated (2MOPP) to all life parts . The applied part is a type B part. The appliance is designed for use in patient environment.

For protection of patients, operators and other persons the following means are implemented:

- An reinforced insulation of input line voltage (life part) from the safety low voltage supply inside the device and the applied part with 2 MOPP
- Basic Insulation (1 MOPP) between safety low voltage and grounded applied part
- Metallic enclosure connected with potential earth, basic insulation (1MOPP) life part to enclosure, potential earth line long enough

For details see also test report EN 60601-1 LGA and insulation diagram QuadroStarPRO.

8.1.3 Protection against electromagnetic phenomena

The following measures were taken:

- EMC-tight metallic enclosure
- internal and external cables are shielded, shields are extensive connected with PE
- line filter for the conductive emission

With these measures the device fulfils the demands of the EN 60601-2:2007 (see also Test Report EN 60601-2 SLG ...).

Form 7301	Risk Management Concept / Report acc. to		 Asclepion Laser Technologies
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 23 of 29 FB 7301 RMR OuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

8.1.4 Protection against laser radiation

For this point see also checklist EN 60601-2-22 (1996) and EN 60825-1.

8.1.5 Protection against unwanted laser emission

The device uses a medical approved foot switch for laser activation. The foot button is mechanical protected against unintentional release, needs a defined mechanical force for activation and has two independent switches implemented.

The device is equipped with a Ready/Standby device.

The device has a laser emergency stop, which disable the laser output immediately.

The device contains of a visible Laser-Ready sign and an optical and acoustical laser emission indicator. The optical indicator is located on the front side of the display.

The laser area can protect by the possibility to connect an external door interlock and/or an external laser warning lamp.

8.1.6 Protection against wrong use

The use of the device is only possible with the key (key switch). The user is informed about the laser danger by the labeling on the device (laser warning label, laser exit aperture, laser emergency button, take notice of user manual) and safety instructions in the user's manual.

8.1.7 Protection against hazardous or wrong output

The control data is calculated by using the treatment parameters set by the user. Variable pulse systems are provided. The drive value is an analog voltage. Its amplitude is proportional to diode current and laser output power. Each laser output (pulsed or continuous) is measured and controlled. If a measured value is lost due to processing steps that are too slow and/or a repetition frequency that is too high, an error message will be generated and the laser operation will be terminated.

The laser output values are controlled by means of the control of the output power of the laser. For this purpose, the laser diode module contains a photo-diode detector. This detector are measured and evaluated. If the difference is less than 20% of the measured value to be expected (stored power adjustment) the laser is in the normal operation mode. If the difference is between 20% and 50% of the measured value to be expected (stored power adjustment) a warning sign will be shown at the display. Otherwise, difference is more than 50%, an error message will be displayed and the laser operation will be stopped.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 24 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

During the commissioning and the final adjustment, the nominal value for the energy control is determined for each device separately and entered into a non-volatile storage of the PEMS.

Adequate instructions for proper maintenance by the user including the demand for the annularly recalibration of the laser output by a well trained service technician ensure the accuracy of the output values during the lifetime of the device.

Master protection:

The measuring system monitors in pulsed and continuous mode permanently derogate energy.

Photodiode over/under energy:

- If the measured power more than 20% from the setpoint, the over/under-energy symbol is displayed. The device is still operational.
- If the measured power more than 50% from the setpoint, the device immediately goes into standby mode. A message over /under energy will be shown.
- A further inspection is carried out in standby mode. System monitors photodiode current in standby mode. If a current in standby measured, then the device is locked with a message. Master set shutdown flag and device shutdown after 5s.

Slave protection:

The slave monitors over energy only. Under energy is not observed in slave.

- If the measured power more than 50% from the setpoint, the device immediately goes to standby mode. A message over /under energy will be send to master. Master show message.
- Slave monitor current in pulsed mode and generate message, if current more than 10% in Ton or Toff time.
- The failure tolerance time of the system is estimate to 500 ms.
- An optical measurement of the laser output
- Laser output energy is strong related to the diode current. Therefore a power adjustment of the output to the diode current is performed. In the used laser parameter field the correlation is sufficient linear. The current is measured by an independent sensor, a shunt located in the diode driver.
- The current signals will checked again the factory calibration curve.
- By controlling the current of the diode and the power output by photodiode we can sure monitor energy deviation of more than 20% of the nominal output.

Form 7301	Risk Management Concept / Report acc. to		 Asclepion Laser Technologies
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 25 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
		Member code/year/month/day-sequential number	

- The most dangerous situation is over-energy by more than 20% of the nominal energy. But also to low energy can result in an insufficient therapeutic effect.

The output energy control is based on the standard EN 60601-2-22:1996 (prEN 60601-2-22:2005) that is interpreted in the following way:

- Section 50 (sec.12.1): Exactness of the operational data
 - a) unit for indicating the value of the laser emission
 - solution/application: 5,7" VGA color display with touch buttons around screen
 - b) indication in SI units
 - solution/application: indication of data: ms / s, Hz, J/cm² / W
- Section 51 (sec.12.4): Protection against hazardous output values
 - a) open-loop or closed-loop admissible
 - solution/application: a open-loop is used for laser.
 - b) Values higher / lower than the intended laser energy are optically or acoustically indicated.
 - solution/application:
 - If the laser energy deviates by more than $\pm 20\%$ in several pulses in succession, an notice sign is shown at the display.
- Section 52 (sec.13): Improper operation and errors
 - a) Output laser energy has been higher than double the set value (+100%) for a period longer than 100ms
 - solution/application: look at c)
 - b) Output of pulsed laser energy, if the laser energy of the previous pulse has exceeded double the set value (+100%).
 - solution/application: look at c)
 - c) Output of repeated pulsed laser energy, if the laser energy of successive pulses has exceeded double the set value (+100%) and the time interval of the successive emissions has been longer than 100ms.
 - solution/application:
 - measurement of the energy values of each pulse
 - measure value generation by peak detection measurement of the photodiode

Solution:

- The measurements are exclusively performed at the ready mode after each pulse with the foot switch pressed.
- The system is checked for:

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 26 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

- o a loss of the measurement data between two pulses
- o an energy deviation higher than 20% (warning sign will be shown)
- o an energy deviation less than 20% (warning sign will be shown)
- o an energy deviation higher than 50% (pulse generation will be stopped immediately and an error information will be shown)
- o an energy deviation less than 50% (pulse generation will be stopped immediately and an error information will be shown)

8.1.8 Failure of mechanical, resistance to heat, expelled parts

The device QuadroStarPRO is designed according to the latest technical standard. All components are fixed. Additional fuses are provided at critical connections (e.g. spring rings or dental disks). All wires and cables are fixed and involved in harness. Sharp corners and edges are protected with edge protection strips. Hoses are secured against slipping out with special circlips. Leakage water can flow away through the holes at the bottom of the housing. The maximum water volume in the system is limited to approximately 200ml. The used hoses and connectors are designed for 6 fold operating pressure. The maximum system pressure is approximately 1 bar. The maximum system pressure (no scanner handpiece connected, self closing water connectors) is limited by a bypass system to 1 Bar.

The device has a metal plate for stability that is inside a closed plastic enclosure. These plastic parts are made from UL listed raw material. The top consists of safety glass with a display behind it.

Heat inside the device

To prevent heat inside the device, it is equipped with acceptable cooling system and automatically control system. The measurements for verify this are carried out by accredited laboratory.

Heat inside the hand piece

The thermoelectric cooler (TEC) for the skin cooling is a source of heat inside the scanner hand piece. This component is driven by high current and can create heat. The total thermal load inside the hand piece is for the TEC maximum 50 Watt. To remove the heat the TEC warm side is water cooled by an air-water heat exchanger. Under normal conditions the cooling capacity is strong enough to prevent overheating. If no scanner handpiece is connected (self closing water connectors) the maximum system pressure is limited by a bypass system to 1 Bar. At this pressure the power consumption of the water pump is less than 5 Watt. Therefore an overheating of the pump is excluded.

In the case of cooler failure a temperature and sensor (monitoring flow) in the water cooling circuit and an additional T-sensor in the scanner hand piece will switch-off the TEC and diode supply immediately.

Form 7301	Risk Management Concept / Report acc. to		
	EN ISO 14971: 2009		
Index: E	Identification number:	MMo/13/10/01-02	Page 27 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc
	Member code/year/month/day-sequential number		

Protection against low temp

The thermoelectric cooler (TEC) for the skin cooling is a source of cooling inside the scanner hand piece. The temperature of the probe is monitored and controlled by a microcontroller system.

Under normal conditions there is a two point regulation of the TEC. In the case of a failure (single fault condition) the surface of the TEC increase above 0°C immediately ($t < 10\text{sec.}$), because of the limited thermal capacity of the cooling area and the cooling power of the TEC.

	Limit	Result
TEC - Nominal "100%"	18°C ± 1K	symmetrical hysteresis 18°C
TEC - Nominal "70%"	22°C ± 1K	symmetrical hysteresis 22°C
TEC - Nominal "30%"	26°C ± 1K	symmetrical hysteresis 26°C
TEC - SkinCooler On (Standby)	$+ 5^\circ\text{C} \leq T \leq \text{TEC - Nominal} + 1^\circ\text{C}$	switch-on the TEC
TEC - Ready Off	$+ 5^\circ\text{C} < T < + 40^\circ\text{C}$	switch-off the TEC and diode supply immediately

Result of risk assessment: For incorrect operation (reduction of temperature) of the TEC occur no hazardous situation and/or no unacceptable risk to patients, operators or third parties. The remaining residual risk is acceptable.

8.1.9 moving parts

There are no moving parts.

8.1.10 Movement

The QuadroStarPRO is a table top device with a total weight under 15kg.

8.1.11 Safety signs

Risks that are not apparent to the operator are indicated by appropriate safety signs. The safety signs are explained in the manual and in the labeling.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 28 of 28 FB 7301 RMR QuPRO V02_MMo131001-02.doc

8.1.12 Cleaning and disinfection, Sterilization and Compatibility

All surfaces of the device are permanently cleanable. Application Parts and parts from the laser handpiece can be disinfected and sterilize during the operational lifetime. The application part is made from stainless steel 1.4301 to reduce allergic reactions.

For more details see Users Manual 149099001 chapter "Cleaning and disinfection, Sterilization" and the cleaning instructions of the handpieces and fibers.

8.1.13 Accuracy and Arrangement of controls and indicators

The device QuadroStar PRO has all indicators and controls in the top. The indicator display is behind a glass panel. The control buttons are around the display. The controls are capacitive touch surfaces. In the top middle of the glass panel is the laser emission LED. All settings must be confirmed. The confirmed settings are released with READY-Button. The indicators are in full decimal places, up to one decimal place. The units are SI units.

8.1.14 Acoustic energy

The device is designed so that human exposure to acoustic energy and vibration cannot result in an unacceptable risk.

The device QuadroStar PRO sends out acoustic energy from the cooling fan. To avoid impair of hearing, the fan and pump noise is reduced to an acceptable level. Check the datasheet of fan and pump.

8.1.15 Batteries

The backup battery of the controller is only available to service engineers with special tools. The discharge current is limited on the circuit board to be 0,3 A. The circuit path is dimensioned to 0,9 A. An additional security is not necessary.

Form 7301	Risk Management Concept / Report acc. to EN ISO 14971: 2009		 Asclepion Laser Technologies
Index: E	Identification number:	MMo/13/10/01-02 Member code/year/month/day-sequential number	Page 29 of 29 FB 7301 RMR QuPRO V02_MMo131001-02.doc

9. Risk from inadequate controlled processes in design, manufacture and maintenance

The processes for design, manufacturing and maintenance are all according to the certified quality system of ALT GmbH acc. EN ISO 13485 and MDD 93/42/EEC Annex II.

Design processes and clinical evaluation were performed by ALT.

Outsourced manufacturing processes (electrical, electronic and mechanical parts) made under the supervision and evaluation of ALT, in accordance to the certified regulations of the quality manual (vendor monitoring).

The manufacturing of electronic PCB's is carried out by external EN ISO 13485 certified manufacturer.

Safety-related power supplies are standard goods of international suppliers, delivered with test certificate, approvals and production control by approved test authorities (such as VDE, KEMA).

Other safety-related parts (power circuit) are a commodity and have the mark of recognized testing laboratories (TUV / VDE / KEMA).

Maintenance and repair processes will be carried out by trained ALT staff or by external contractors trained and certified by the ALT service.

These processes are all designed in a way, that they can be classified as appropriate as part of this risk management process and the remaining residual risk is acceptable.

10. Risk control and risk verification

The verification and validation of the risk measures was carried out according to the process instructions PA 735 "Design Verification" and PA 736. "Design validation".

FB 7351e	Verification / Validation Report	
Revision Index: B	Identification number: ZI/13/03/04-01 Member Code/year/month/day-sequential no.	Page 1 of 3 FB 7351e UsabilityValid V01 ZI130304-01.doc

Verification	☐	Validation	☒ Usability
Product / part:	QuadroStarPRO	Material number:	1490ff
Cause:	☒ Product development	☐ Product change	☐
Original document:	☒ Requirement specification of: FB 7321 LH QuPRO V02 ZI121122_01 in KA 5.1	☐ Design specification of:	☐
	☐ Risk analysis of:	☐ Change-application no.:	☐

Serial number: *see test reports*

Brief description: (traceability of testing)

This validation is based on the forerunner device QuadroStar, which is broadly similar to the QuadroStarPRO for the user. Besides the new technique inside the device (which is not relevant for the user) and the new shape the differences for the user are the yellow wavelength as a new option and the new display.

The intended use, parameters, handling and accessories didn't change to the QuadroStar version so far, which never caused problems in recent years and was used with great success.

The validation test can therefore be focused to the new wavelength and the new display.

Selected measuring and test equipment: n/a

Documents of the enclosure: (test reports etc.)

Protocol of assessment of case series of yellow wavelength by Dr. Wood: TP 7342 Prot Wood 310113

Protocols of user evaluation: ZI/13/04/09

Final assessment of results: (brief assessment, summary measures)

The first test treatments of the optional available new yellow wavelength showed the advantage compared to green as expected.

As required is the operation of the device easy to understand and intuitively.

☐ VERIFICATION / VALIDATION NEGATIVE ☒ VERIFICATION / VALIDATION POSITIVE

if negative, derived measures:

Plan Ident-No. ZI/13/02/28/-01 drawn by:

14.1.13 Zimmermann *Zimmermann*
Date Name Signature

Verification/Validation performed:

9.4.13 Zimmermann *Zimmermann*
Date Name Signature

Plan checked and released:

14.1.13 Katzer *Katzer*
Date Name Signature

Verification/Validation checked and released:

9.4.13 Katzer *Katzer*
Date Name Signature

Cause of latest modification FB: review					
FB drawn by process owner R&D:			FB checked and released by QA:		
05.07.2012	I. Streit	sdg. Streit	06.07.2012	R. Thieme	sdg. Thieme
Date	Name	Signature	Date	Name	Signature

FB 7351e	Verification / Validation Report	
Revision Index: B	Identification number: ZI/13/03/04-01 Member Code/year/month/day-sequential no.	Page 2 of 3 FB 7351e UsabilityValid V01 ZI130304-01.doc

History:

Rev	Brief description of change	Date	Name	Signature
01	First edition			

Linked documents:

document	abbreviation	file name
Protocol of assessment of case series of yellow wavelength by Dr. Wood		TP 7342 Prot Wood 310113
Protocols of user evaluation		TP 7342 Prot UserValid_Vs1_Zi 130409-1

NOTE:

The intention for the development of the new QuadroStar version – QuadroStarPRO – was mainly to implement the technical improvements of the last years and to reduce the production costs. The intended use including operator profile etc., principle of action to the tissue, physical functionality and the handling / POFs are identical to the QuadroStar. Therefore the long experience with the QuadroStar and forerunners of more than 20 years in our company is the basis for the intended improvements.

As an example: the new yellow wavelength was not possible until now with respect to the compactness and weight of the intended device and costs also. It is identical in all aspects mentioned above. The advantage is a reduction of necessary Fluence for the treatment and therefore a certain reduction of the normal side effects due to the slightly higher blood absorption and slightly lower melanin absorption compared to green (known from other, large and expensive yellow lasers including forerunner devices of ALT).

All of the possible hazards which could be caused by the device in question are identical to our medical device QuadroStar, which to our best knowledge never caused health problems beside the normal side effects stated in the user manual (documented also in App3-BfArM to ClinEval QuadPRO_ZI_120829-01 in KA 7.6). Therefore these hazards have to be managed for the new device in the same way as for the forerunner device.

The novel for the user is besides the remodeled, more trendy shape the display and its handling only. Of course the specification for this is based on our experience with all our medical laser devices.

Therefore the validation of the usability is focused on the difference: yellow wavelength and display

Validation

No.	Description	Expected result	Result	Acceptance Yes(v) No(-)
A	Validation of yellow wavelength			
1.	Fluence needed for successful treatment (compared to green)	Less Fluence	About 40% less Fluence	yes
2.	Less side effects compared to green	Less side effects	Less pain Less side effects	yes
3.			Yellow more effective	yes
B	Validation of new display	(validation by school grades 1 very good - 5 bad)		
4.	Switching on uncomplicated	Better than 3	1.1	yes
5.	Display bright enough		1.6	Yes
6.	All details can be seen clearly, also from side		1.6	Yes
7.	Display to see with safety goggle also		2.1	Yes
8.	Clear separation of display / setting		1.4	Yes

FB 7351e	Verification / Validation Report	 Asclepion Laser Technologies
Revision Index: B	Identification number: ZI/13/03/04-01 Member Code/year/month/day-sequential no.	Page 3 of 3 FB 7351e UsabilityValid V01 ZI130304-01.doc

No.	Description	Expected result	Result	Acceptance Yes(v) No(-)
9.	Clarity	Better than 3	1.7	Yes
10.	All important infos are shown		1.8	Yes
11.	Icons intuitively to understand		2.3	Yes
12.	Operation structure logically, self-explaining		2.7	Yes
13.	Operation fast and safe		2.6	Yes
14.	Operation guidance intuitively		2.6	Yes
15.	Short learning, easy to remember		2.7	Yes
16.	High safety of operation, min. mistakes		2.6	Yes
17.	Setting of parameters		1.6	Yes
18.	Buttons are easy to operate		3	Yes
19.	Warning sign, alert good to recognize		1.4	Yes
C	Evaluation of manual			
20.	Manual good to understand		1.2	Yes
21.	Warnings/ hints clear and precise		1.2	Yes
22.	Good description of operation		1.2	Yes
23.	Good description of cleaning and maintenance		1.2	Yes

Assessment

- The advantage of treatments with yellow light compared to green could be shown very clear as expected on the basis of literature and own experience (see clinical evaluation in KA 7.7)
- The new display was accepted by this group of 7 persons who were not involved in the development of the device. The evaluation was better than 3 = middle, except no. 18: operation of buttons – this will be corrected by adjustment of sensitivity of the buttons.
- The manual was evaluated by the group as very good.

CL 7351	RL 93/42/EWG * Conformity Evaluation Essential Requirements acc. to CD 93/42/EEC *	Asclepion Laser Technologies
Revision index: C	Identification number: MMo/14/02/20-01 Identification code/year/month/day-number	page 1 of 32 CL 7351 Conf. MDD 93_42_EEC QuPRO V02 MMo130220-01.doc

Tested device: **QuadroStarPRO**
Name and/or model: **1490-1495**
Seriennummer:
Serial number: *n./a.*
Herstellungsdatum (wenn bekannt): *n./a.*
Date of manufacturing (if known):
Hersteller
Manufacturer: **Asclepion Laser Technologies GmbH**
Anschrift
Address: **Brüsseler Straße 10, 07747 Jena, Deutschland**
Applied standard: Richtlinie 93/42/EWG vom 14.06.1993 über Medizinprodukte *
 Council Directive 93/42/EEC of 14 June 1993 concerning medical devices *
Attachments:
Note:
Examined by: **M. Möller**
 Name Date *21.02.2014*
Checked by: **I. Streit**
 Name Date *21.02.2014*

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Revision index: C

Identification number: MMo/14/02/20-01

Identification code/year/month/day-number

page 2 of 32
CL 7351 Conf. MDD 93_42_EEC QuPRO
V02 MMo130220-01.doc**Identification****Details of examiner:**

Identification of the person responsible for examining the product under inspection:

Name: MICHAEL MÖLLER
print full name (block letters)Position: Project manager R&D
print full title

If the above named person is not an employee of the manufacturer of the product to be examined, state the details of the examiner's employer or organization:

Organisation: _____

Address: _____


Medizinprodukt/Medical device

Details of product to be inspected:


User manual

1490 99 001_Vs2-10.Febr14-GA_QuPRO-DE

 Asclepion Laser Technologies GmbH 07747 Jena, Brüsseler Straße 10, Germany	
Model / SN / 00100	QuadroStarPRO 100-240 V ~ 50/60 Hz 100 V ~ 635 VA 240 V ~ 780 VA
 CE 1275  15°C 	

Picture of the device

Scanner handpiece

RL 93/42/EWG *

Conformlty Evaluation Essential Requirements acc. to CD 93/42/EEC *

Asclepion
Laser Technologies

Revision index: C

Identification number: MMo/14/02/20-01

Identification code/year/month/day-number

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CL 7351

RL 93/42/EWG *

Conformity Evaluation Essential Requirements acc. to CD 93/42/EEC *

Asclepion
Laser Technologies

Revision index: C

Identification number: MMo/14/02/20-01

Identification code/year/month/day-number

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Nr.	zutreffend	Richtlinie/ Norm/ Regeln	Zugehörige CL	Vorhanden
30	☐	UL 60601-1 (U.S. national deviations) <i>Medical Electrical Equipment, Part 1</i> <i>General Requirements for Safety</i>		
31	☐	21 CFR 1040 (04-1-02 edition) / Notice 50 FDA <i>Code of federal regulations - Food and Drugs 21.</i>		
32	☐	CAN/CSA-ISO 13485:2003 <i>Quality management systems- Requirements for regulatory purpose</i> CMDCAS <i>Canadian Medical Devices Conformity Assessment System</i>		
33	☐	ISO 14971: 2007 <i>Medical devices - Application of risk management to medical devices</i> <i>Medical devices - Application of risk management to medical devices.</i>		
35	☒	DIN EN 207 <i>Persönlicher Augenschutz.</i> <i>Personal eye-protection equipment.</i>		
36	☒	DIN EN 50 419: 2006-06 <i>Kennzeichnung von Elektro- und Elektronikgeräten.</i> <i>Marking of electrical and electronic equipment.</i>		
37	☒	EN ISO 14971: 2009 <i>Medizinprodukte: Anwendung des Risikomanagements auf Medizinprodukte.</i> <i>Medical devices - Application of risk management to medical devices.</i>		

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 1 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

Verification	☐	Validation	☒ Literature route
Product / part:	QuadroStarPRO	Clinical Evaluation	
		Material number:	149x
Cause:	☐ Product development	☒ Product change	☐
Original document:	☒ Requirement specification of: FB.7321 LH Taurus V01 ZI/12/03/20-01	☒ Design specification of: FB7322 Pflichtenheft QuadPRO V01 MMo120626-01	☐
	☐ Risk analysis of:	☐ Change-application no.:	☐
Serial number	N/A		

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The goal of the new version of the QuadroStar⁺ was to improve the technology and to reduce the production costs - including the possibility to implement this reliable technology for the yellow wavelength, which has the same intended use as the green wavelength and is more selective (but so far more expensive). Therefore there was no change relating to the intended use as described in the User Manual of the QuadroStar so far. This was the request of the Requirement Specification (Lastenheft ZI/12/03/20-01) and was realized in the development of the new version. There is no need for a clinical study.

The wavelengths of the device in question are standard laser treatments since more than 30 years with own forerunner devices since more than 25 years. This assessment includes the evaluation of literature of lasers with the same specs like the laser QuadroStarPRO.

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
Based on the evaluation of the information arising by the analysis of the literature the final evaluation, that the benefit of the new device and the application is much higher than the risk for the patient and the user. The expected side effects and risks are marginal and not higher than the risks of other devices (including forerunner devices of the QuadroStarPRO), which are used since years for identical indications.

☐ VERIFICATION / VALIDATION NEGATIVE ☒ VERIFICATION / VALIDATION POSITIVE
 If negative, what are the resulting measures:

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
Reason for latest change in form sheet: Initial Compile Process					
Form sheet compiled by process owner R&D:			Form sheet reviewed & approved by QM:		
17.07.05	I. Streit	Sgd. Streit	10.09.05	R. Thieme	Sgd. Thieme
Date	Name	Signature	Date	Name	Signature

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technology
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 2 of 19 FB 7351 KlinBew QuadroStarPRO ver. 1.0 ZI120822-01.doc

History:

Rev	Brief description of change	Date	Name	Signature
01	First edition	4.12.12	G.Zimmermann	<i>G. Zimmermann</i>

Clinical Evaluation

- 0 Preliminary note
- 1 Qualification of the author and the reviewer
- 2 Argument for literature route
- 3 Evaluation of literature
 - 3.1 Equivalence of the devices used in the studies
 - 3.2 Study design and ranking
 - 3.3 Evaluation of the scientific studies regarding effectiveness
 - 3.3.1 Effectiveness demonstrated for yellow wavelength
 - 3.3.2 Effectiveness demonstrated for green wavelength
 - 3.3.3 Effectiveness demonstrated for NIR wavelengths
 - 3.4 Evaluation of the scientific studies regarding side effects
 - 3.4.1 Side effects reported for yellow wavelength
 - 3.4.2 Side effects reported for green wavelength
 - 3.4.3 Side effects reported for NIR wavelengths
 - 3.5 Evaluation of contraindications
 - 3.5.1 Contraindications for yellow and green wavelengths
 - 3.5.2 Contraindications for NIR wavelengths
 - 3.6 Evaluation of customer experience
 - 3.7 Evaluation of reports of the BfArM
- 4 Final assessment of results

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 3 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

Preliminary note

The goal of the development of a new version of the QuadroStar⁺ – the QuadroStarPRO - is the technological improvement by keeping the intended use, which is proved to be very successful. The technology is now more than 10 years old and the components available today allow to reduce the size of the device and production costs. Additionally there is now the possibility to implement this reliable technology for the yellow wavelength, which has the same intended use as the green wavelength and is even a little more selective (but so far more expensive – therefore the less expensive green wavelength will be kept as an option).

As with the QuadroStar⁺ so far the QuadroStarPRO will have different device versions, which are identical besides the wavelength. The green and infrared wavelengths of the new device are the same as the QuadroStar so far.

The device in question will have a yellow wavelength additionally, which we already had with similar devices in the past (copper bromide lasers ProYellow and YellowStar) for the identical intended use of the green, but stopped selling because of technical reasons. With the new technological developments it is possible to add a version with a yellow wavelength and the same design and reliability.

All wavelengths of the device in question have been used for laser treatments since more than 30 years now with own forerunner devices since about 25 years and are standard in the meantime.

Because of higher blood and melanin absorption the yellow and the green wavelength are more selective in view of treatments of vascular and pigmented lesions compared to the near infrared (NIR) wavelengths of the device and therefore mainly used for precise skin treatments.

The NIR on the other hand has a deeper penetration because of less scattering and is therefore used for a wide range of applications in small surgery mainly, comparable to electrosurgery but more precise.

Because of the relatively long pulse duration of milliseconds up to seconds or even continuous emission (cw) the tissue is coagulated and with higher intensities and/or absorption vaporized.

Please see also:

App. 1 of Customer Requirement Specification: Specification of Application (ZI/12/03/20-01)

Table 1: Laser devices which are used for the intended use of the QuadroStar

A. For green and yellow wavelengths,
intended use: superficial vascular, benign pigmented lesions and other skin irregularities

Type of device	Wavelength	comment
Argon lasers Own devices: DLS 5/DL 5000, 1983 - 1998	514nm (green) (+ blue 488nm)	Are used in dermatology since more than 30 years Are replaced now by newer technologies Relatively high risk of side effects due to low power and therefore longer pulse duration
Dye lasers	577 – 590nm (yellow)	Are used in dermatology since more than 25 years advantage: yellow, which is a bit more selective than green Disadvantage: Very large devices, high running costs, therefore mainly used in hospitals or specialized laser centers
Copper bromide Copper vapour lasers Own devices: ProYellow, YellowStar 1998 - 2002	578nm (yellow) (+ green 511nm)	Are used in dermatology since more than 15 years advantage: yellow, which is a bit more selective than green Disadvantage: large devices, technology is not very reliable

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 4 of 19 FB 7351 KinBew QuadroStar ZI120822-01.doc

Type of device	Wavelength	comment
KTP / LBO laser Own devices: BeautyStar QuadroStar, QuadroStar+ Since 2002	532nm (green)	Are used in dermatology since more than 10 years advantage: small, reliable, table top devices

B. For NIR wavelengths,
intended use: coagulation and vaporization of tissue

Type of device	Wavelength	comment
Nd:YAG lasers Own devices: JenLas MED 100, JenLas MED 60, Ergolas 1991 - 1998	1064nm	Are used in small surgery since more than 30 years Are replaced now by newer technologies (diode lasers)
Diode lasers Own devices: QuadroStar, QuadroStar+ Since 2004	800 – 980nm	Are used in small surgery since more than 15 years advantage: small, table top devices

1. Qualification of the author and reviewer

The author of this evaluation, Mrs. G. Zimmermann is the Senior Product Manager and the person in duty for medical validation and application by Asclepion Laser Technologies with a lot of experience to carry out clinical evaluations.

Mrs. G. Zimmermann has a diploma in physics and has experience in application of medical lasers for more than 20 years, especially for surgical lasers, lasers for dermatology and aesthetics including the QuadroStar and devices using intense pulsed light (IPL), radiofrequency and acoustic waves.

After graduation she developed laser devices as a scientific engineer for more than 10 years. Since 1988 she was responsible for medical application of different lasers including QuadroStar and worked together with physicians in various specialties. She also performed laser safety courses including medical application for physicians for 18 years now and educated the sales staff and doctors in using the devices.

She was in charge to supervising own application studies and evaluating scientific literature. Working as Product Manager for more than 20 years she was also responsible for fixing of application parameters of new lasers as well as supporting the development in respect to ergonomically questions, writing the application part of user manuals, evaluating new devices in respect to medical use and ergonomics.

The reviewer of this evaluation, Mr. André Steps, is not an employee of the manufacturer Asclepion Laser Technologies GmbH.

Mr. Steps has a prediploma in medicine and a sport science degree. Since more than 3 years he educates doctors and nurses in the theoretical principles and practical application of lasers, IPLs and other devices for dermatological and esthetic treatments. In this time he did a vast number of trainings all over the world.

Therefore he has an excellent overview about all relevant lasers and devices used worldwide for the intended use of the QuadroStar as well.

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 5 of 19 FB 7351e KinBew QuadPRO V01 ZI120822-01.doc

2. Argument for literature route

The goal of the new version of the QuadroStar was to improve the hard- and software of the device, which is some years old, and to reduce the production costs.

The effectivity of treatments with the QuadroStar and forerunners (see table 1) is proven, also compared to similar devices. Equally important, during all the time there were no complaints about serious damages to health. The number of side effects (small blisters, crusts and temporary pigment changes), which cannot be excluded totally, was very low.

Additionally the remaining risks are identical to the current QuadroStar, which did not cause any harm to patients or users in the past years besides in some cases the possible side effects stated in the manual.


Therefore the intended use (including indications and target group/contraindications) of the new version QuadroStarPRO is requested to be identical compared to the current version of the QuadroStar (Requirement Specification ZI/10/11/18-01). The application parameters related to the clinical performance, benefit/risk ratio, safety and risks are the same: the range of wavelength, of fluence - respectively power - and of pulse duration as well as the scanner and the skin cooling system.

The only difference is the additional yellow wavelength, but:

- the intended use including all details mentioned above is identical to the green
- besides the wavelength all application parameters are the same as for the green
- we already had in the past devices with this yellow wavelength, so we have own positive experience (ProYellow, YellowStar)
- the yellow is also well established in the market since more than 25 years and even a little better compared to the green (because more selective, which causes less side effects):

Fig. 1:

Blood absorption
for green and yellow
wavelengths

GREEN 532nm:

very high blood absorption, very low water absorption (not shown)
therefore very selective specially for vascular treatments

YELLOW 577nm:

at maximum of blood absorption, most selective absorption
even a little better for vascular treatments

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Page 6 of 19 FB 7351 KlinBew QuadPRO ZI120822-01.doc
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	

Due to the little differences between the green and yellow wavelengths, the intended use is identical including indication (superficial vascular, benign pigmented lesions and other skin irregularities), contraindications, target group, restrictions to user, application technique etc.

Also the different wavelengths of the NIR have nearly the same tissue effect of coagulation and vaporization and are therefore also identical in view of aspects mentioned above.

The clinical data were gained from scientific studies, textbooks and clinical guidelines. Since the indicated treatments are established techniques (textbooks and guidelines available) and the present device is a further development of a long history product line of Asclepion Laser Technologies (ALT), sufficient clinical data with appropriate evidence is available and renders a clinical study unnecessary.

The assessment by literature route is possible.

3. Evaluation of literature

The skin treatments which are performed with the yellow and green wavelength, as well as the indications treated with the NIR are standard laser treatments since a long time. Therefore there are clinical guidelines (App.4 and 5 for yellow/green lasers) and textbooks (e.g. [1] for NIR lasers) documenting that these kind of treatments are effective, safe and to the best advantage for the patients. This literature is used besides the own experience with forerunner devices as a basis for this clinical evaluation and for the application part of the User Manual. For details of indications, contraindications and side effects stated in this literature with comments to utilization in the User Manual see App. 7.

The objective of this clinical evaluation is to show the effectivity and safety of treatments with the device in question using actual literature, guidelines and textbooks. A literature search was performed to find actual published valid data. The overall results of the literature search of studies documented in App.2 are listed below. These articles will serve as basic information for the clinical evaluation of the QuadroStarPRO in general while the outcomes for the single wavelengths will be discussed separately.

The overall articles provide clinical, histological as well as metadata findings. All the selected clinical data is assessed concerning its informative value using an evidence-based ranking system recommended by the Agency for Healthcare Research and Quality. This scale is most suitable to appraise the quality of interventional procedures and hence makes it possible to organize the studies due to their degree of quality on the basis of a widely recognized norm (Kunz R, Ollenschläger G, Raspe H, Jonitz G, Donner-Banzhoff, 2007: Lehrbuch Evidenzbasierte Medizin in Klinik und Praxis. Dt Ärzte-Verlag, Köln).

As a result of this rating procedure, most of the interventional studies feature the evidence class IIa and IIb. Lower evidence levels were only accepted if the article presented possible adverse events. The particular classification is stated in the table below. No further gradation is undertaken.

The information content varies from study to study. Essential data, such as description of the used laser device, quantitative analysis of the outcomes, declaration of exclusion criteria and side effects is stated adequately whereby the outcome of each study can be related internally to its treatment procedure allowing a reliable performance- and risk analysis. The study populations and respective treatment outcomes are summarized for similar indications. The overall study population covers all Fitzpatrick skin types I-VI.

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FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 7 of 19 FB 7351 KlinBew QuadroPRO V01 ZI120822-01.doc

3.1 Comparability of the devices used in the studies

All the articles were carefully screened to be in line with the intended clinical application of the QuadroStarPRO built on the experience gathered with the QuadroStar⁺. With respect to the technical accordance of the devices used in the studies, it can be said that all studies apply devices equivalent to the QuadroStarPRO, which are already listed in table 1. In the guidelines (App.4 and 5) as well as in the textbook [1] these devices are summarized to be equivalent for the treatments in question. Therefore there is no need to evaluate the parameters of all the devices applied in the studies in detail. The listed datasets can therefore be considered as suitable for the demonstration of performance and safety.

3.2 Study design and ranking

Table 2 Studies for Clinical Evaluation of QuadroStarPRO^{YELLOW}

Author	Title	Year	Journal/Proceedings	Evidenc class
Studies supporting performance by presenting histological findings				
(Neumann, et al., 1992)	Comparative histochemistry of port-wine stains after copper vapor laser (578 nm) and argon laser treatment.	1992	J Invest Dermatol. Vol. 99(2), pp. 160-167	Ila
(Chung, et al., 1996)	Histological responses of port wine stains in brown skin after 578 nm copper vapor laser treatment	1996	Lasers Surg. Med. Vol. 18(4), pp. 358-366	Iib
(Neumann, et al., 1993)	Results and tissue healing after copper-vapour laser (at 578 nm) treatment of port wine stains and facial telangiectasias.	1993	Br J Dermatol. Vol. 128(3), pp. 306-312	Iib
(Walker, et al., 1989)	Histology of port wine stains after copper vapour laser treatment	1989	British Journal of Dermatology Vol. 121(2), pp. 217-223	Iib
Studies supporting performance by presenting clinical findings				
(Waner, et al., 1993)	A comparison of copper vapor and flashlamp pumped dye lasers in the treatment of facial telangiectasia.	1993	J Dermatol Surg Oncol. Vol. 19(11), pp. 992-998	Ila
(Chung, et al., 1997)	Copper vapour laser treatment of port-wine stains in brown skin	1997	Australasian Journal of Dermatology Vol. 38(1), pp. 15-21	Iib
Karsai, S., Czarnecka, A. & Raulin, C.	Treatment of Xanthelasma Palpebrarum Using a Pulsed Dye Laser: A Prospective Clinical Trial in 38 Cases	2010	Dermatologic Surgery Vol. 36(5), pp. 610-617	Iib
(Kaufman, et al., 1995)	Treatment of adenoma sebaceum with the copper vapor laser	1995	Journal of the American Academy of Dermatology Vol. 33(5, Part 1), pp. 770-774	Iib
(Key J, 1992)	Selective destruction of facial telangiectasia using a copper vapor laser	1992	Archives of Otolaryngology-Head & Neck Surgery Vol. 118(5), pp. 509-513	Iib
(Madan, et al., 2009)	Pulsed-dye laser treatment of angioma serpiginosum	2009	Clinical and Experimental Dermatology Vol. 34(5), pp. 186-188	Iib
(McCoy, et al., 1996)	An evaluation of the copper-bromide laser for treating telangiectasia.	1996	Dermatol Surg. Vol. 22(6), pp. 551-557	Iib
(McCoy, 1997)	Copper bromide laser treatment of facial telangiectasia: Results of patients treated over five years	1997	Lasers Surg. Med. Vol. 21(4), pp. 329-340	Iib

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 8 of 19 FB 7351 KinBew QuadPRO Z1120822-01.doc

(Park, 2008)	Pulsed dye laser treatment for viral warts: A study of 120 patients	2008	The Journal of Dermatology Vol. 35(8), pp. 491-498	IIb
(Park, et al., 2007)	Pulsed Dye Laser Therapy for Pediatric Warts	2007	Pediatric Dermatology Vol. 24(2), pp. 177-181	IIb
(Sadick, 2002)	The Utilization of a New Yellow Light Laser (578 nm) for the Treatment of Class I Red Telangiectasia of the Lower Extremities	2002	Dermatologic Surgery Vol. 28(1), pp. 21-25	IIb
(Sheehan-Dare, 1994)	Copper vapour laser (578 nm) and flashlamp-pumped pulsed tunable dye laser (585 nm) treatment of port wine stains: results of a comparative study using test sites.	1994	Br J Dermatol. Vol. 130(4), pp. 478-482	IIb
(Pickerling, et al., 1990)	Copper vapour laser treatment of port-wine stains and other vascular malformations	1990	British Journal of Plastic Surgery Vol. 43(3), pp. 273-282	III
(Tay, 2012)	Treatment of infantile hemangiomas with the 595-nm pulsed dye laser using different pulse widths in an Asian population.	2012	Lasers Surg Med. Vol. 44(2), pp. 93-66	III
Studies suitable for appraising side effects				
(Sheehan-Dare, 1994)	Copper vapour laser (578 nm) and flashlamp-pumped pulsed tunable dye laser (585 nm) treatment of port wine stains: results of a comparative study using test sites.	1994	Br J Dermatol. Vol. 130(4), pp. 478-482	IIa
(Karsai, et al., 2010)	Treatment of Xanthelasma Palpebrarum Using a Pulsed Dye Laser: A Prospective Clinical Trial in 38 Cases	2010	Dermatologic Surgery Vol. 36(5), pp. 610-617	IIb
(Kaufman, et al., 1995)	Treatment of adenoma sebaceum with the copper vapor laser	1995	Journal of the American Academy of Dermatology Vol. 33(5, Part 1), pp. 770-774	IIb
(Key J, 1992)	Selective destruction of facial telangiectasia using a copper vapor laser	1992	Archives of Otolaryngology-Head & Neck Surgery Vol. 118(5), pp. 509-513	IIb
(Madan, et al., 2009)	Pulsed-dye laser treatment of angioma serpiginosum	2009	Clinical and Experimental Dermatology Vol. 34(5), pp. 186-188	IIb
(McCoy, et al., 1996)	An evaluation of the copper-bromide laser for treating telangiectasia.	1996	Dermatol Surg. Vol. 22(6), pp. 551-557	IIb
(McCoy, 1997)	Copper bromide laser treatment of facial telangiectasia: Results of patients treated over five years	1997	Lasers Surg. Med. Vol. 21(4), pp. 329-340	IIb
(Neumann, et al., 1993)	Results and tissue healing after copper-vapour laser (at 578 nm) treatment of port wine stains and facial telangiectasias.	1993	Br J Dermatol. Vol. 128(3), pp. 306-312	IIb
(Park, 2008)	Pulsed dye laser treatment for viral warts: A study of 120 patients	2008	The Journal of Dermatology Vol. 35(8), pp. 491-498	IIb
(Park, et al., 2007)	Pulsed Dye Laser Therapy for Pediatric Warts	2007	Pediatric Dermatology Vol. 24(2), pp. 177-181	IIb
(Sadick, 2002)	The Utilization of a New Yellow Light Laser (578 nm) for the Treatment of Class I Red Telangiectasia of the Lower Extremities	2002	Dermatologic Surgery Vol. 28(1), pp. 21-25	IIb

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 9 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

(Pickering, et al., 1990)	Copper vapour laser treatment of port-wine stains and other vascular malformations	1990	British Journal of Plastic Surgery Vol. 43(3), pp. 273-282	III
(Adamic, et al., 2007)	Vascular lasers and IPLS: Guidelines for care from the European Society for Laser Dermatology (ESLD)	2007	J Cosmet Laser Ther Vol. 9(2) Journal of Cosmetic and Laser Therapy, pp. 113-124	IV

Table 3 Studies for Clinical Evaluation of QuadroStarPRO^{GREEN}

Author	Title	Year	Journal/Proceedings	Evidence class
Studies supporting performance by presenting clinical findings				
(Erceg, et al., 2010)	A Comparative Study of Pulsed 532-nm Potassium Titanyl Phosphate Laser and Electrocoagulation in the Treatment of Spider Nevi	2010	Dermatologic Surgery Vol. 36(5), pp. 630-635	Ib
(Collyer, et al., 2010)	Comparison of Treatment of Cherry Angiomata With Pulsed-Dye Laser, Potassium Titanyl Phosphate Laser, and Electrodesiccation A Randomized Controlled Trial	2010	Archives of Dermatology Vol. 146(1), pp. 33-37	Ila
(Kundu, et al., 2009)	Comparison of Electrodesiccation and Potassium-Titanyl-Phosphate Laser for Treatment of Dermatos Papulosa Nigra	2009	Dermatologic Surgery Vol. 35(7), pp. 1079-1083	Ila
(ÖÖzden, et al., 2011)	Clinical comparison of potassium-titanyl-phosphate (KTP) versus neodymium:YAG (Nd:YAG) laser treatment for lower extremity telangiectases	2011	Journal of Dermatological Treatment Vol. 22(3), pp. 162-166	Ila
(Uebelhoer, et al., 2007)	A Split-Face Comparison Study of Pulsed 532-nm KTP Laser and 595-nm Pulsed Dye Laser in the Treatment of Facial Telangiectasias and Diffuse Telangiectatic Facial Erythema	2007	Dermatologic Surgery Vol. 33(4), pp. 441-448	Ila
(Grunewald, et al., 2008)	Treatment of superficial vascular skin changes with a 532 nm lithium-borate laser	2008	Medical Laser Application Vol. 23(2), pp. 55-60	Ilb
Studies suitable for appraising side effects				
(Erceg, et al., 2010)	A Comparative Study of Pulsed 532-nm Potassium Titanyl Phosphate Laser and Electrocoagulation in the Treatment of Spider Nevi	2010	Dermatologic Surgery Vol. 36(5), pp. 630-635	Ib
(Kundu, et al., 2009)	Comparison of Electrodesiccation and Potassium-Titanyl-Phosphate Laser for Treatment of Dermatos Papulosa Nigra	2009	Dermatologic Surgery Vol. 35(7), pp. 1079-1083	Ila
(Uebelhoer, et al., 2007)	A Split-Face Comparison Study of Pulsed 532-nm KTP Laser and 595-nm Pulsed Dye Laser in the Treatment of Facial Telangiectasias and Diffuse Telangiectatic Facial Erythema	2007	Dermatologic Surgery Vol. 33(4), pp. 441-448	Ila
(Grunewald, et al., 2008)	Treatment of superficial vascular skin changes with a 532 nm lithium-borate laser	2008	Medical Laser Application Vol. 23(2), pp. 55-60	Ilb
(Dawson, et al., 2007)	Adverse Events Associated With Nonablative Cutaneous Laser, Radiofrequency, and Light-Based Devices	2007	Seminars in Cutaneous Medicine and Surgery Vol. 26(1), pp. 15-21	IV

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 10 of 19 FB 7351 KInBew QuadPRO V01 ZI120822-01.doc

Table 4 Studies for Clinical Evaluation of QuadroStarPRO^{810,940,980}

Author	Title	Year	Journal/Proceedings	Evidence class
Studies supporting performance by presenting clinical findings				
(Leclère, et al., 2010)	A prospective randomized study of 980 nm diode laser-assisted venous ulcer healing on 34 patients	2010	Wound Repair and Regeneration Vol. 18(6), pp. 580-585	Ib
(Nordon, et al., 2011)	A prospective double-blind randomized controlled trial of radiofrequency versus laser treatment of the great saphenous vein in patients with varicose veins.	2011	Ann Surg. Vol. 254(6), pp. 876-881	Ib
(Crisan, et al., 2008)	Use of laser-therapy for treatment of oral and maxillofacial hemangioma.	2008	Rev Med Chir Soc Med Nat Iasi. Vol. 112(2), pp. 512-516	Ila
(Dziecluchowicz, et al., 2011)	A Prospective Comparison of Four Methods of Endovenous Thermal Ablation	2011	Polish Journal of Surgery Vol. 83(11), pp. 597-605	Ila
(Desmytère, et al., 2010)	Endovenous Laser Ablation (980 nm) of the Small Saphenous Vein in a Series of 147 Limbs with a 3-Year Follow-up	2010	European Journal of Vascular and Endovascular Surgery Vol. 39(1), pp. 99-103	Ilb
(Desmytère, et al., 2007)	Endovenous 980-nm laser treatment of saphenous veins in a series of 500 patients	2007	Journal of Vascular Surgery Vol. 46(6), pp. 1242-1247	Ilb
(Gibson, et al., 2007)	Endovenous laser treatment of the short saphenous vein: Efficacy and complications	2007	Journal of Vascular Surgery Vol. 45(4), pp. 795-803	Ilb
(Park, et al., 2008)	Endovenous Laser Ablation of the Incompetent Small Saphenous Vein with a 980-nm Diode Laser: Our Experience with 3Years Follow-up	2008	Eur J Vasc Endovasc Surg Vol. 36(6), pp. 738-742	Ilb
(Trelles, et al., 2012)	Laser-assisted lipolysis in the treatment of gynecomastia: a prospective study in 28 patients	2012	Lasers in Medical Science, pp. 1-8	Ilb
(Vuylsteke, et al., 2008)	Endovenous Laser Treatment of Saphenous Vein Reflux: How Much Energy Do We Need to Prevent Recanalizations?	2008	Vascular and Endovascular Surgery Vol. 42(2), pp. 141-149	Ilb
Studies suitable for appraising side effects				
(Nordon, et al., 2011)	A prospective double-blind randomized controlled trial of radiofrequency versus laser treatment of the great saphenous vein in patients with varicose veins.	2011	Ann Surg. Vol. 254(6), pp. 876-881	Ib
(Dziecluchowicz, et al., 2011)	A Prospective Comparison of Four Methods of Endovenous Thermal Ablation	2011	Polish Journal of Surgery Vol. 83(11), pp. 597-605	Ila
(Desmytère, et al., 2010)	Endovenous Laser Ablation (980 nm) of the Small Saphenous Vein in a Series of 147 Limbs with a 3-Year Follow-up	2010	European Journal of Vascular and Endovascular Surgery Vol. 39(1), pp. 99-103	Ilb
(Desmytère, et al., 2007)	Endovenous 980-nm laser treatment of saphenous veins in a series of 500 patients	2007	Journal of Vascular Surgery Vol. 46(6), pp. 1242-1247	Ilb

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 11 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

(Gibson, et al., 2007)	Endovenous laser treatment of the short saphenous vein: Efficacy and complications	2007	Journal of Vascular Surgery Vol. 45(4), pp. 795-803	IIb
(Park, et al., 2008)	Endovenous Laser Ablation of the Incompetent Small Saphenous Vein with a 980-nm Diode Laser: Our Experience with 3Years Follow-up	2008	Eur J Vasc Endovasc Surg Vol. 36(6), pp. 738-742	IIb
(Trelles, et al., 2012)	Laser-assisted lipolysis in the treatment of gynecomastia: a prospective study in 28 patients	2012	Lasers in Medical Science, pp. 1-8	IIb
(Vuytstele, et al., 2008)	Endovenous Laser Treatment of Saphenous Vein Reflux: How Much Energy Do We Need to Prevent Recanalizations?	2008	Vascular and Endovascular Surgery Vol. 42(2), pp. 141-149	IIb
(Theivacumar, 2009)	Arterio-Venous Fistula Following Endovenous Laser Ablation for Varicose Veins	2009	Eur J Vasc Endovasc Surg Vol. 38(2), pp. 234-236	III
(Yildirim, et al., 2009)	Treatment of an Unusual Complication of Endovenous Laser Therapy: Multiple Small Arteriovenous Fistulas Causing Complete Recanalization	2009	CardioVascular and Interventional Radiology Vol. 32(1), pp. 166-168	IV

3.3 Evaluation of the scientific studies regarding effectiveness

3.3.1 Effectiveness demonstrated for YELLOW wavelength

The studies reflect the major application area of the yellow wavelength laser achieving convincing results in the treatment of vascular lesions and other skin irregularities. All studies dealing with treatment of port wine stains and telangiectasia showed good to excellent results in the majority of the treated subjects.

Particularly the study of Sue McCoy (1997) is important for this evaluation, because it reports about the huge number of 570 patients treated with our own forerunner device: 87.4% had good to excellent result, only 12.6% had less than 50% clearance of telangiectasias.

Histologic examinations revealed that the thermal injury zone was limited to blood vessels and perivascular collagen layers leaving overlying connective tissue as well as hair follicles and glands viable (Neumann, et al., 1992); (Neumann, et al., 1993). So it clearly demonstrates in practice the selectivity of the wavelength used, which was seen in theory.

Clinical data suggests that the effectiveness of a vascular treatment with yellow laser light mainly depends on the vessel color and diameter as well as the skins melanin level (Chung, et al., 1997); (Sadick, 2002).

Histologic findings demonstrated that the additional absorption in melanin leads to less thermal injury in darker skin even with fluences around 10 J/cm² (Chung, et al., 1996) causing a higher rate of poor treatment outcomes (Chung, et al., 1997).

Regarding the blood vessel characteristics it can be stated that the treatment of superficial dark vessels with a diameter up to 1 mm resulted in higher success rates than other deeper, brighter or larger vessels (Pickering, et al., 1990; Neumann, et al., 1993). Larger vessels require higher fluences and/or pulse widths for clearance which may result in a higher risk for complications (Sadick, 2002).

When compared to alternative therapy options, the copper bromide technology achieved significantly lower fading scores than pulsed dye laser technique in port wine stains (Sheehan-Dare, 1994).

Park et al. demonstrated the possibility to treat non-pigmented lesions such as warts with a total clearance rate of about 50 %.

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 12 of 19 FB 7351 KinBew QuadroPRO V01 ZI120822-01.doc

The following treatment outcomes were achieved in the selected articles representing an overall population of about 1300 treated subjects:

Table 5 Studies demonstrating performance of yellow wavelength

Citation	Indication	Results
(Chung, et al., 1996)	port wine stains	morphological changes in basal cells appeared starting at 6 J/cm ² ; loss of viability of the epidermal cells starting at 8 J/cm ² ; in dermis, damage confined to blood vessels and surrounding collagen; at 10 J/cm ² a wedge-shaped diffuse coagulation necrosis; vascular selectivity without epidermal damage cannot be achieved at 578 nm in the brown skin of Koreans
(Chung, et al., 1997)	port wine stains	n = 18 (38.2%) with good to excellent results; n = 13 (27.6%) with fair results, n = 16 (34.0%) with poor results; doses above 10 J/cm ² with wedge-shaped diffuse coagulation necrosis and a marked loss of viability of epidermis and underlying dermis
(Karsai, et al., 2010)	xanthelesma palpebrarum	approximately 2/3 of the lesions with clearance greater than 50%; 1/4 with clearance greater than 75%
(Kaufman, et al., 1995)	adenoma sebaceum	good to excellent cosmetic results in 100% of treated patients
(Key J, 1992)	telangiectasia	n = 18 with satisfactory results; n = 2 (patients with telangiectasia associated with acne rosacea) with unsatisfactory response
(Madan, et al., 2009)	angioma serpiginosum	after Ø3.75 treatments (range 1–7) complete clearance in five lesions, excellent response (> 80% improvement) in four, and good response (50–80%) in three lesions
(McCoy, et al., 1996)	telangiectasia vascular Neavi	n = 1 with excellent 90-100% clearance with 2 treatments; n = 9 with 75-90% clearance with Ø2.1 treatments; n = 5 with 50-75% clearance; n = 1 with 25-50% clearance with 1 treatment; n = 1 with <25% clearance with 1 treatment
(McCoy, 1997)	telangiectasia	n = 72 (12.6%) with poor result (0-50% clearance) with Ø1.5 treatments; n = 99 (17.4%) with good result (50-75% clearance) with Ø1.85 treatments; n = 399 (70%) with excellent result (>75% clearance) with Ø2.17 treatments; cheek vessels with higher clearing rate than nasal alae and/or tip vessels; exposure time and energy density necessary increased with vessel diameter
(Neumann, et al., 1992)	telangiectasia	histology 15 min post: with up to 12 J/cm ² epidermal cells viable; thermal injury confined to blood vessels in papillary and mid-reticular dermis; perivascular collagen damaged; overlying connective tissue viable; hair follicles, sebaceous and sweat glands not affected; vessel part facing the laser beam damaged; at least 20 J/cm ² : full thickness necrosis; depth of injury correlated with pulse duration and energy density; not significant between 8-12 J/cm ²
(Neumann, et al., 1993)	telangiectasia port wine stains	port wine stains: n = 6 (24%) with excellent, n = 7 (28%) good, n = 7 (28%) fair results; n = 5 (20%) showed no response; with sub-threshold doses for clinical whitening (vessel-selective damage) only moderate blanching in one pink and one red PWS; other lesions with higher, non-specific energy densities; best results in light-red to red PWS; purple colour with fair or no response; facial telangiectasias: n = 15 (45%) with excellent results, n = 8 (24%) with good results, n = 4 (12%) with fair results, n = 6 (19%) with no improvement; up to 12 J/cm ² damage was confined to the vessel walls and a cuff of perivascular collagen; >15 J/cm ² diffuse arc-shaped, non-specific epidermal and dermal coagulation necrosis, sharply demarcated from viable surrounding tissue
(Park, 2008)	viral warts	49.5% overall clearance rate: flat warts (67.6%), peri-ungual warts (51.1%), plantar warts (47.6%) and common warts (44.3%)
(Park, et al., 2007)	pediatric warts	48.1% clearance rate (99/206 treated warts eradicated); effective treatment response (≥50% clearance) in 60.7% of treated patients; effective fluences were between 8.5 and 9.5 J/cm ²
(Pickering, et al., 1990)	port wine stains	n = 36 completed treatment with removal of ≥ 70% of lesion; n = 10 received ≥ 4 four treatments with ≥ 30% reduction in area and slight/no colour change; n = 190 are continuing treatment (outcome between good and poor); most of colour/area reduction occurs with first few treatments
(Sadick, 2002)	telangiectasia	75-100% clearance in 71.8% of patients after Ø1.7 treatments; 50-75% clearance in 13% of patients; 25-50% clearance in 8.8% of patients; <25% clearance in 4.3% of patients; in 2.1% of patients worse or no change; n = 2 without significant lightening; 76.1% of patients displayed great satisfaction
(Sheehan-Dare, 1994)	port wine stains	the PTDL produced significantly better fading scores as determined by clinical assessment and reflectance spectrophotometry compared with the CVL
(Tay, 2012)	infantile hemangioma	short and long pulse durations equally effective; shorter pulse durations with slightly higher incidence of side effects

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 13 of 19 FB 7351_KlinBew_QuadPRO V01 ZI120822-01.doc

3.3.2 Effectiveness demonstrated for GREEN wavelength

As the absorption coefficient of the green laser wavelength in hemoglobin is only slightly lower than the yellow one, the statements mentioned in the subchapter above can be adopted for the QuadroStarPRO^{GREEN} also.

Collyer et al and Erceg et al. showed the ability of the green wavelength laser to achieve outcomes equal to standard surgical procedures while causing less pain. Analogue to the yellow laser wavelength, the vessel size should not exceed a certain diameter. Öözden recommends to treat vessels larger than 1 mm with laser systems like the Nd:YAG (Öözden, et al., 2011).

Table 6 Studies demonstrating performance of green wavelength

Citation	Indication	Results
(Collyer, et al., 2010)	cherry angioma	in pairwise comparisons, cherry angiomas treated with electrodesiccation were significantly less improved than were those receiving KTP laser (P=.003) and with PDL (P=.001); the effects of KTP laser and PDL on texture were not different (P=.50).
(Erceg, et al., 2010)	spider nevi	no statistically significant difference between electrocoagulation and KTP-laser treatment; tendency toward clinical outcomes in favor of the KTP; KTP laser treatment was less painful
(Grunewald, et al., 2008)	telangiectasia; spider naevi; hemangioma; xanthelasma; adenoma sebaceum; spider veins	successful treatment of telangiectasia with very high patient satisfaction; spider naevi cured within 4 weeks with 2 treatments; immediate clearing of senile hemangioma in most regions; no success on spider veins even with highest energy and pulse width; n= 2 with xanthelasma without clearing; good clearance of adenoma sebaceum in one patient after various laser treatments
(Kundu, et al., 2009)	papulosa nigra	no treatment difference for subjective effectiveness (p = .06) or subjective confidence improvement between electrodesiccation and KTP-Laser; significant treatment difference for subjective discomfort (p = .002) in favor of KTP
(Öözden, et al., 2011)	telangiectasia	after the third treatment session, average improvement scores of 2.44, 1.31 and 3.75, 3.23 were given for the KTP and Nd:YAG laser-treated sides, respectively; the KTP laser has very low efficacy with vessels larger than 1 mm and should not be elected when treating such vessels

3.3.3 Effectiveness demonstrated for NIR wavelengths

The studies below demonstrate the successful treatment of 1445 patients with various indications, mainly endovascular laser treatment of varices (EVLV). The latter is a rather new example for a NIR laser treatment in the field of so called 'small surgery'.

The very broad range of indications in the different specialties is treated by precise coagulation or vaporization of tissue. As already mentioned above, these treatments are standard since more than 30 years and described in detail in textbooks for the different specialties.

An overall standard book for laser treatments in different specialties is quoted therefore in the application part of the User Manual:

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Table 7 Studies demonstrating performance of NIR wavelengths

Citation	Indication	Results
(Crisan, et al., 2008)	hemangioma, low-flow vascular malformations	volume regression of lesions and reduction of vascular signal as measured ultrasonographically varied between 45% and 95%
(Desmytère, et al., 2007)	EVLV GSV	successful occlusion, defined as vein occlusion with absence of flow, was noted in 501 GSVs (98%) in the immediate postoperative period; this initial success rate remained almost constant during the follow-up period to reach 97.1% at 4 years; at 1 year, a complete disappearance of the GSV or minimal residual fibrous cord was noted for all patients

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 14 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

(Desmytère, et al., 2010)	EVLt SSV	the initial technical success rate was 100%; SSV remained closed in 114 of 117 limbs (97%) after 1 year, all of 61 limbs after 2 years and all of 30 limbs after 3 years; three SSVs with recanalization
(Dzieciuchowicz, et al., 2011)	EVLt	technical success in 99% of limbs; the differences between the groups did not reach statistical significance; recanalization rates were 7% after EVLA980, 2% after EVLA810 and 0% after RFA and EVLA1470
(Gibson, et al., 2007)	EVLt GSV, SSV	procedure was technically successful in 100% of limbs; at 2-4 d post, 12 (5.7 %) limbs were found to have a deep venous thrombosis (DVT); no DVT was present in any limb at final follow-up scan 2 to 11 months post
(Leclère, et al., 2010)	venous ulcers	at the ninth follow-up visit, in the control group, the ulcers had decreased on an average to 94.3% of the original area, in the laser group, the decrease was to 74.2%; this difference was not significant
(Nordon, et al., 2011)	EVLt GSV	at 3 months, occlusion was 97% for RFA and 96% for EVLT; after EVLT, the median time for return to work was 7 days (1-60), compared to 9 days (0-28) after RFA; neither modality proved to be superior; RFA is associated with less procedural pain, analgesic requirement, and bruising
(Park, et al., 2008)	EVLt SSV	the technical success rate was 100%; the SSV remained closed in 89 of 93 limbs (96%) 1 month post, all of 82 limbs after 6 months, 77 limbs after 1 year, 71 limbs after 2 years and 55 limbs after 3 years; 4 limbs with recanalization
(Theivacumar, 2009)	EVLt GSV	All 3 patients showed the expected reduction in varicose vein score; concomitant venous and arterial wall thermal injury or needle trauma during administration of tumescent anaesthesia may cause this rare complication
(Trelles, et al., 2012)	lipolysis	at 6 months, 18 patients (64.3%) scored the result as "very good", 6 as "good" (21.4%), 3 as "fair" (10.7%) and 1 as "poor" (3.6%); physicians scored the photographs as 'very good' in 22 patients (78.6%), "good" in 5 (17.9%), and "fair" in 1 (3.6%)
(Vuylsteke, et al., 2008)	EVLt GSV	the clinical success rate was 100%; at 6 months, recanalization occurred decreasing occlusion rate to 90.6%; authors found 52 J/cm2 mean fluence in the occlusion group and 43.7 J/cm2 in the non-occlusion group; the occlusion rate on long term is fluence dependent
(Yildirim, et al., 2009)	EVLt GSV	2nd endovenous treatment; recanalization persisted; surgical treatment; 60 J/cm did not decrease vessel diameter due to the large diameter of GSV; AV fistula formation may be related to number of thrombi produced during ablation; suggest higher energy doses (>80 J/cm)

3.4 Evaluation of the scientific studies regarding side effects

Due to the long pulses (up to continuous) of this type of laser the emission of this device causes thermal effects – coagulation up to vaporization with higher intensities, which is similar for all wavelengths of the device in question. Therefore the side effects which could be caused in general are also similar.

The difference is the penetration depth, which is nearly the same for yellow and green, but larger for NIR. Also the higher melanin absorption of yellow and green compared to NIR as well as the type of treatment (skin treatments vs. small surgical treatments mainly) result in some differences.

3.4.1 Side effects reported for YELLOW wavelength

Particularly the study of Sue McCoy (1997) is important for this evaluation, because it reports about the huge number of 570 patients treated with our own forerunner device: only 3 patients had temporary hypopigmentation, no lasting adverse effects were seen.

The overall risk profile illustrates safe treatment characteristics within the limits of the recommended use as stated in the user manual.

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 15 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

It seems adequate to consider Fitzpatrick skin type IV to VI as contra-indications as evidence for a higher risk of depigmentation exists.

This was stated in the manual as well as tanned skin and photosensitivity (also by medication).

The side effects reported in the studies of the table below are known and stated in the manual* (subchapter of application 'possible side effects'):

- Immediately after treatment:
erythema (redness) and edema (swelling) in the treatment area for a few hours or days, but can resist up to weeks
- Sometimes blisters, crusts and temporary pigment changes
- Relatively seldom: small sunken scars
- In single cases: formation of Granyloma pyogenicum in the treated area
- Little risk: superinfections

* Note: Purpura and bruising is not mentioned in the manual, because it is a side effect regularly seen using short pulsed dye lasers only, because related to their much shorter pulse duration (about 100 times shorter than that of the device in question), that's **why** more mechanical effect: disruption of small vessels instead of coagulation by heat as it is with all the other lasers including the device in question.

All detailed description of side effects in the manual are based on the following standard literature:

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Details:

The studies suitable for the evaluation of side effects represent an overall of about 1200 exposed patients covering all skin types I-VI. The follow up duration was at least 4 weeks. The most common laser treatment related side effect is transient erythema and swelling which must be seen as a normal part of the reorganisation process of the skins microstructure. For purpura and bruising see 'Note' above.

Pigment changes were observed in about 10% of the exposed subjects: While on the one hand post-inflammatory hypo- or hyperpigmentation could be the result of overtreatment, the outcomes of Chung et al. and Tay et al. indicate a direct relation between patient's skin type and risk of pigment changes. However, all pigment changes subsided within at most six month and no persistent hyper- or hypo-pigmentation was reported in any case.

Furthermore, the effectiveness of this wavelength in relation to the acceptable risks seems to be limited to vessel diameters of up to 2 mm (Sadick, 2002).

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 16 of 19 FB 7351 KInBew QuadPRO V01 ZI120822-01.doc

Table 8 Summary of yellow laser related side effects

Citation	Side effects	Follow up [months]
(Adamic, et al., 2007)	pain; purpura; bruising; bleeding; hematoma; disruption; swelling; discoloration; blisters; scabs; infection; reactivation of herpes simplex on the face; hyperpigmentation; hypopigmentation; skin texture changes; scarring; lesion persistence	n.a.
(Chung, et al., 1997)	temporary slight hyperpigmentation in nearly all patients (resolved within 6-12 weeks)	3
(Karsal, et al., 2010)	purpura, swelling in 100% of lesions (for 10-12 d, 2-5 d); hyperpigmentation in 3 lesions	1
(Kaufman, et al., 1995)	superficial crusting 2-3 d post; n= 1 with transient hypopigmentation (resolved within 4 weeks);	n.a.
(Key J., 1992)	edema (subsided within 3-4 d); about 50% with scabbing at day 3 post until day 12; n= 3 (15%) with post-inflammatory hyperpigmentation (treatment with hydroquinone; resolved within 6-8 weeks)	4
(Madan, et al., 2009)	n= 8 with purpura (subsided in 7-10 d); n= 2 with transient hyperpigmentation (resolved after 6 months)	n.a.
(McCoy, et al., 1996)	erythema (for 24 h); slight edema (for up to 3 d); crusting	2
(McCoy, 1997)	n= 3 with hypopigmentation	2
(Neumann, et al., 1993)	reddening; edema; n= 13 with transient hyperpigmentation (resolved within 6-12 weeks); n= 18 with slight atrophic scarring	2
(Park, 2008)	minimal-to-mild pain during treatments in most patients; minimal to mild post-treatment discomfort; crusting (resolved in 2 weeks); 7 drop out due to intolerable pain at 10 J/cm ²	3-36
(Park, et al., 2007)	crusting (resolved in 2 weeks); 1 drop out due to intolerable pain during treatment; n= 9 with moderate to severe pain during treatment, posttreatment discomfort, minimal scarring, hypopigmentary changes	1
(Pickering, et al., 1990)	n= 6 with scarring (less than 5% of treated area was affected); n= 2 with hyperpigmentation; n= 2 with hypopigmentation	n.a.
(Sadick, 2002)	n= 2 (4.4%) with mild epidermal crusting; n= 3 (6.5%) with temporary hyperpigmentation	3
(Sheehan-Dare, 1994)	n= 4 with slight atrophy; n= 1 with moderate atrophy; 4 sites with hyperpigmentation	4
(Tay, 2012)	n= 2 with blistering; n=3 with hyperpigmentation; n= 4 with hypopigmentation; n= 3 with textural changes; edema;	n.a.

3.4.2 Side effects reported for GREEN wavelength

The general aspects are the same as for the yellow wavelength. The difference is minor due to the little difference in blood/melanin absorption.

The overall risk profile illustrates safe treatment characteristics within the limits of the recommended use as stated in the user manual.

Table 9 Summary of green laser related side effects

Citation	Side effects	Follow up [months]
(Dawson, et al., 2007)	transient erythema; edema; purpura; transient-linear crusting; temporary dyspigmentation; rarely scarring or textural changes	1
(Erceg, et al., 2010)	pain was reported as only adverse event	3
(Grunewald, et al., 2008)	erythema in 80% of cases for 24 h; no pigment changes; crusting above 20 J/cm ²	1.5
(Kundu, et al., 2009)	n= 14 with pain, redness; n= 1 with scarring; n= 2 with swelling; 1 drop out after the first treatment because of hyperpigmentation, which resolved completely within 3 months	-

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 17 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

3.4.3 Side effects reported for NIR wavelengths

NIR lasers of the QuadroStarPRO type are mainly used for small surgical and endoscopic treatments. Coagulation and vaporization of tissue lead to side effects like pain, inflammation and swelling in and around the treated area like with other methods (e.g. electrosurgery – but no side effects related to the electric field of these devices).

These side effects and the coagulated tissue respectively the lack of tissue after vaporisation have to be balanced versus the advantages of the laser treatment compared to conventional surgery (no bleeding, short hospitalization etc.) in each individual case by the physician. This is explicitly stated in the User Manual.

As an example the side effects of endovascular treatments are stated below.

Table 10 Summary of NIR-laser related side effects

Citation	Side effects	Follow up [months]
(Desmytère, et al., 2007)	no DVT; ecchymoses in 60% for Ø 2 weeks; transitory paresthesia in 7% of treated legs for Ø 2 weeks; moderate pain in 9.3%	6
(Dzieduchowicz, et al., 2011)	bruising (in 27% of patients); tightness or pain (13%); paresthesia (4.2%)	36
(Gibson, et al., 2007)	n= 6 with temporary paresthesia or hypoesthesia (mostly related to the phlebectomies); n= 6 with residual skin pigmentation; n=2 with active venous ulcer (healing at 1 month post)	6
(Leclère, et al., 2010)	discomfort when moving the arms; slight darkening from bruising	6
(Park, et al., 2008)	no deep venous thrombosis; ecchymoses in 60% with a median duration of 2 weeks; temporary paraesthesia (mostly hypoaesthesia) observed in 40% of treated legs with a median duration of 2 weeks	36
(Theivacumar, 2009)	In one patient from EVLA980 group the procedure was not performed because of laser generator failure; overall incidence of deep venous thrombosis 1.3% (2/154) and increased to 4.9% (2/41) in 810 nm group; n= 12 with saphenous nerve paresthesia on the 10th postoperative day	12
(Trelles, et al., 2012)	n= 3 (<0.2%) with arterio-venous fistula; 2 fistulae closed spontaneously with conservative management	3-18
(Vuytsteke, et al., 2008)	multiple small arteriovenous fistula; nearly complete recanalization	-
(Yildirim, et al., 2009)	postoperative pain was considerably less after RFA, 2.6% with skin burn; 1.3% with paresthesia 2.6% with thrombophlebitis	3

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	 Asclepion Laser Technologies
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 18 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

3.5 Evaluation of contraindications

3.5.1 Contraindications for yellow and green wavelengths

The contraindications stated in the User Manual are based on our own experience with forerunner devices and double-checked with the actual guidelines for care for these treatments of the German and European Society for Laser dermatology.

Contraindication User Manual QuPRO	Contraindication European Guideline	Notes
Pregnancy	No impact on pregnancy Laser treatment should be discontinued in the first 2 trimesters May be postponed until after delivery	No clinical studies
Breastfeeding		No clinical studies
Children		Almost all treatments can be postponed urgent treatments under responsibility of the physician
Tanned skin	Patients with a tan should delay treatment	
Skin type IV-VI	For darker skin more care	Because of high melanin absorption high risk of side effects
Light sensitivity		Could lead to more side effects
Taking of medication which increases light sensitivity		Comparable to tan and dark skin
Cancer. Precancerosis, NCN in the treatment area		Often discussed at congresses
Kasabach-Merrit-Syndrom	Kasabach-Merrit-Syndrom	new adopted as warning!!
Herpes simples in the treatment area		Treatment can trigger Herpes, also mentioned in guideline
Tendency of hypertrophic scarring or formation of keloids	Individuals with history of hypertrophic scarring or keloids should be treated with caution	

Additionally it is stated that it is the responsibility of the physician to decide the individual case.

3.5.2 Contraindications for infrared wavelengths

As already stated in the chapter about side effects it is the responsibility of the physician to decide the individual case, therefore it is also stated in the User Manual.

3.6 Evaluation of customer experience

The direct forerunner device is sold since about 10 years now, the green and yellow wavelengths even much longer.

For effectiveness and safety see the studies of Sue McCoy with a forerunner device mentioned in 3.3.1 and 3.4.1.

In all these years we didn't had any complaint about treatments with this device beside the (minor) side effects stated in the User Manual.

3.7 Evaluation of reports of the BfArM

As shown in App. 3 there were no reports or recalls relevant for the device in question.

FB 7351e	Verification / Validation Report Clinical Evaluation QuadroStarPRO (name of project: TAURUS)	
Revision Index: A	Identification Number: ZI/12/08/22-1 Shortname/Year/Month/Day-Ser. No.	Page 19 of 19 FB 7351 KlinBew QuadPRO V01 ZI120822-01.doc

4. Final assessment of results

The indications stated in the Manual are standard laser treatments since a long time and published for the skin treatments (yellow and green wavelength) in the European and German guidelines (App. 4 and 5) respectively for small surgery e.g. in the chapter 'Therapeutische Leitlinien zur Laserbehandlung' of the textbook [1].

The existence of these guidelines and standard textbooks for different specialties proves already in its own that the laser treatments in question are to the best advantage for the patients balancing the effect versus the minor side effects.

Additional the evaluation of studies realized with forerunner devices of our own company or other devices which are equivalent showed the efficacy of the treatments indicated in the intended use of this new QuadroStarPRO laser device.

Side effects of the skin treatments due to individual overtreatment can be blisters, crusts, pigment changes and scars, which are small, because limited to the target (e.g. vessel). They are normally transient despite the small scars, which are seldom. Side effects cannot be totally avoided, but can be reduced to a minimum by doing a test treatment as described in the user manual.

The contraindications found in the literature are included in the manual as shown in detail in the related chapter of this evaluation.

The user is informed about remaining risks in the user manual:

- the possible side effects in chapters 6.3.7 and 6.4.4
- the contraindications in chapter 1.3.2

The indications, contraindications and side effects are checked for completeness - see in App. 6: Application CL Manual QuPRO.

All warnings and informations necessary in the User Manual due to the Risk analysis of the new device are checked for completeness by App. 8: Risk control by Manual QuPRO

The conformity with all requirements of the Medical Device Directive is checked by App.7: CL 7356 Checkliste Gebrauchsanweisung QuadroStarPRO.

Final assessment:

Based on this analysis the final evaluation is, that the benefit of the new device and the application is much higher than the risk for the patient and the user. The expected side effects and risks are marginal and not higher than the risks of other devices (including the older versions of the QuadroStarPRO) which are used since years for the indications in question.

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