Cordis Corporation

14201 NW 60th Avenue, Miami Lakes, FL 33014

Mailing Address: P.O. Box 025700, Miami, FL 33102

T (786) 313 2000 F (786) 313 2080

April 20, 2017

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Cordis diagnostic catheters are available in a broad variety of French sizes and configurations. These catheters combine an atraumatic tip with either a braided or non-braided body.

Cordis diagnostic catheters are designed to deliver radiopaque contrast medium to selected sites in the vascular system.

Unknown.

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Possible complications include, but are not limited to the following:

- air embolism
- hematoma at the puncture site
- infection
- perforation of the vessel wall.

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Introduce the catheter into the blood vessel using a vascular entry technique of choice. PTFE coated guidewires are recommended for use with Cordis diagnostic (angiographic) catheters.

Specifications

The following specifications are identified on the catheter hub (refer to Figure 1):

1. maximum pressure rating (psi)
2. catheter outer diameter (F)
3. catheter length (cm)
4. maximum guidewire diameter (in.)

The color of the hub printing indicates the catheter French size as follows:

- 4F (1.35 mm) Red
- 5F (1.65 mm) Grey
- 6F (2.0 mm) Green
- 7F (2.3 mm) Orange
- 8F (2.7 mm) Blue

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein. Descriptions or specifications in Cordis printed matter, including this publication, are

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meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

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Cordis diagnostic catheters are available in a broad variety of French sizes and configurations.

Cordis Catheter Extensions are available in different configurations.

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Cordis Catheter Extensions can be used to connect the angiographic catheter to the power injector. Cordis Catheter Extensions are available in different configurations.

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When necessary to extend catheter please attach the end of the extension catheter with the metal spinning device to the proximal end of the catheter (hub) by turning the spinner clockwise until it stops. To detach the extension catheter rotate the spinning device counterclockwise.

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The manufacturer warrants conformance of this medical device to the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer's use instructions are met.

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein. Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties.

Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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Cordis medical devices should be stored at temperatures ranging from +15°C to + 25°C, in a place protected from direct sunlight.

Cordis medical devices being in operation are stable in temperature conditions from +32 °C to +42 °C and to biological liquids.

Sterility maintaining term - 3 years.

After use all used components and packing materials may be potentially biologically dangerous. Their treatment and disposal are performed according to the accepted medical practice, and also according to local, state and federal laws and provisions. The products which had a contact with blood and / or body fluids are subject to disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class). Unused products (which had no contact with blood and / or body fluids), including those with expired shelf-life, are subject to disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

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Instruction for Use Diagnostic Catheters

Cordis Angiographic Catheters with Marker Bands are available in 5 French. The Cordis Angiographic Catheters with Marker Bands contain multiple radiopaque markers, circumferentially mounted on the catheter. For number of markers refer to product label.

Cordis Angiographic Catheters with Marker Bands are designed to provide angiographic visualization and linear measurement of the vasculature when combined with the delivery of radiopaque contrast media to selected sites in the vascular system.

Unknown.

Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure. Possible complications include, but are not limited to, the following: air embolism, hematoma at the puncture site, infection, thrombosis, hemorrhage, dissection, perforation or other damage of the vessel wall. Movement of the marker bands along the catheter can result in inaccurate reference and device sizing. Dislodgement of the marker bands into the vascular system can result in additional intervention, embolism, thrombosis or other vascular complications.

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ic heparinization.

7. Avoid excessive tension on the device during manipulation. Extreme care to avoid stretching or elongation must be exercised during manipulation and withdrawal.
8. If resistance is felt during manipulation, determine the cause of resistance before proceeding and confirm SUPER TORQUE MB Angiographic Catheter positioning under high quality fluoroscopic observation.
9. Measurement: The distance "X" in mm between two subsequent marker bands is in accordance with the value on the label, the distance should be read as shown in Figure 2.

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may come off the catheter and dislodge into the vascular system.

- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- Do not expose to organic solvents.
- Do not exceed maximum pressure rating printed on label and hub.

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- Avoid entrapment of the catheter between other endovascular devices and the vessel wall.
- Avoid excessive friction on the catheter; avoid simultaneous introduction of the catheter and aortic graft devices through the same sheath.
- Store in cool, dark, dry place.
- Do not use if the package is open or damaged.
- Do not use the catheter if the "Use By" date on the package label has expired.
- Do not resterilize.
- Exposure to temperatures above +54°C may damage the catheter.

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The manufacturer warrants conformance of

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this medical device with the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer's use instructions are met.

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Cordis medical devices should be stored at

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Instruction for Use Guiding Catheters

The Cordis Guiding Catheter has a large nontapered lumen that allows contrast medium injections and pressure monitoring and facilitates the intravascular passage of interventional devices.

The guiding catheter is intended for use for intravascular introduction of interventional/diagnostic devices into the coronary or peripheral vascular systems.

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- air embolism
- hematoma at the puncture site
- infection
- perforation of the heart
- vessel damage, dissection or perforation
- vasospasm
- ischemia
- hemorrhage
- arrhythmia
- reaction to contrast media
- death

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1. Remove the guiding catheter from its packaging.

2. Inspect the guiding catheter upon removal from packaging to verify that it is undamaged.

Warning: Do not use a guiding catheter that has been damaged in any way. If damage is detected, replace with an undamaged guiding catheter.

3. Prior to use, flush the guiding catheter lumen with a heparinized saline solution.

4. Introduce the guiding catheter into the vasculature through the catheter sheath introducer, and/or over an indwelling guidewire using a percutaneous entry technique of choice.

5. Set up a continuous heparinized saline flush through the sidearm of a hemostatic valve attached to the guiding catheter hub. Note: It is recommended that a continuous heparinized saline flush be maintained between the guiding catheter and any intraluminal device passed through it.

6. Under fluoroscopic guidance, advance the guiding catheter over the guidewire or introducer until desired position is attained.

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or to the introduction of other intravascular devices or infusion of contrast agents.

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Risk of reuse: This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use. Reuse of this product, including after reprocessing and/or resterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

Do not use with Ethiodol or Lipiodol (Ethiodol and Lipiodol are trademarks of Guerbet S.A.) contrast media, or other such contrast media which incorporates the components of these agents.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein. Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the prod-

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uct at the time of manufacture and do not constitute any express warranties.

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Introducer Guide

The Introducer Guide has a large nontapered lumen that allows for intravascular passage of catheters and infusion of heparinized normal saline. The Introducer Guide, with sideport tubing and valve, has a pre-shaped distal tip section with radiopaque BRITE TIP. The Introducer Guide is pre-shaped to facilitate its positioning. Package with Introducer Guide includes vessel dilator. A vessel dilator is provided to facilitate the percutaneous entry of the Introducer Guide by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

The Introducer Guide is intended for use for intravascular introduction of interventional and/or diagnostic devices into the peripheral vasculature.

None known.

Procedures requiring percutaneous catheter and/or sheath introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

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- air embolism
- hematoma at the puncture site
- infection
- intima tear/dissection
- perforation of the vessel wall
- thrombus formation.

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place.

5. Follow standard practice for vessel puncture. Insert a guidewire of appropriate length and diameter.
6. Using a gentle twisting motion, advance the Introducer Guide/vessel dilator assembly over the guidewire and through the tissue into the vessel. Continue to advance the assembly until the distal tip is in the desired position.
7. Holding the Introducer Guide in place, withdraw the dilator and, if applicable, the guidewire.
8. Follow standard practice of using a continuous drip of anticoagulant fluid under pressure through the sideport when the Introducer Guide is in the vessel.
9. Follow manufacturer's recommendations for the catheter or device being introduced via the Introducer Guide. Caution: Verify on the label that the Introducer Guide internal diameter (ID) is of proper size to ensure compatibility with the inserted device.
10. When removing the Introducer Guide, re-insert the dilator and wire fully into the Introducer Guide to aid in straightening the tip portion. Thereafter, remove the dilator, wire, and Introducer Guide as a unit.

Caution: Upon withdrawal of dilator, catheter, or other inserted device, aspiration through the sideport is recommended to remove any fibrin accumulation.

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Discard the Introducer Guide after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Introducer Guides are extremely difficult to clean after exposure to biological materials and may cause adverse patient events if reused. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from reuse of the Introducer Guide.

Do not use with Ethiodol or Lipiodol (Ethiodol and Lipiodol are trademarks of Guerbet S.A.) contrast medium, or other such contrast media that incorporates the components of these agents.

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Cordis medical devices should be stored at temperatures ranging from +15°C to +25°C, in a place protected from direct sunlight.

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
These medical devices are for single use only and are not subject to maintenance and repair.

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Authorized Representative of the manufacturer in the Russian Federation:
 Johnson & Johnson LLC
 Russian Federation, Krylatskaya street 17,
 building 2, Moscow, 121614
 Tel.: + 7 (495) 580 77 77;
 Fax: + 7 (495) 580 78 78

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Instruction for Use Guiding Catheter

The Cordis OPTEASE Retrieval Catheter is an 80 cm long, 10F catheter, with radiopaque tip.

The OPTEASE Retrieval Catheter is intended for the percutaneous removal of the Cordis OPTEASE

Retrievable Vena Cava Filter from the inferior vena cava (IVC).

The OPTEASE Retrievable Filter can be retrieved within a specified period after implantation (refer to the OPTEASE Retrievable Filter Instructions for Use) or remain implanted as a permanent filter.


The self-centering OPTEASE Retrievable Filter is laser cut from nickel titanium alloy (Nitinol) tubing. The proximal and distal baskets of the OPTEASE Retrievable Filter consist of struts in a six diamond shape configuration. A single row of fixation barbs is present at the cranial end of the struts. These barbs, intended for fixation to the vessel wall, are extensions of the parallel struts. A hook is centrally located at the caudal basket extremity and allows for filter retrieval using a snare (Figure 1).

Retrieval of the OPTEASE Retrievable Filter is from femoral venous access only.

Recommended accessories to the OPTEASE Retrieval Catheter are listed in Table 1.

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
1. Microvena Snare (wire loop) correctly engaged to the filter retrieval hook. These are in line with the OPTEASE Retrieval Catheter. Maintain a distance of approximately 5 mm between the tip of the Microvena Snare Catheter and the OPTEASE Retrieval Filter retrieval hook

2. Microvena Snare Catheter with radio-paque marker
3. OPTEASE Retrieval Catheter with radi-opaque tip

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1. Microvena Snare (wire loop) incorrectly engaged to the OPTEASE Retrieval Filter retrieval hook.
2. Microvena Snare Catheter with radio-paque marker
3. OPTEASE Retrieval Catheter with radi-opaque tip.

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Advancing the OPTEASE Retrieval Catheter to collapse the OPTEASE Retrieval Filter.

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Table 1 (EN).

Recommended Accessories to the OPTEASE Retrieval Catheter

Recommended Accessories

Dimension/Description

Cordis J-tip EMERALD 0.035" (0.89 mm) Guidewire

260 cm length

Cordis AVANTI+ Catheter Sheath Introducer

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The Cordis OPTEASE Retrieval Catheter is intended for the percutaneous removal of the Cordis OPTEASE Retrieval Vena Cava Filter from the inferior vena cava (IVC)..

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- Retrieval of the filter for patients with an on-going high risk for pulmonary embolism (PE).
- Retrieval of the filter with significant amounts of trapped thrombus within the filter that present a risk for embolism during retrieval.

Procedures requiring percutaneous interventional techniques should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during the procedure.

Possible retrieval procedure complications include, but are not limited to, the following:

- air embolism
- hematoma at the puncture site
- rupture of the vessel wall
- perforation of the vessel wall
- intimal tear
- inability to retrieve the filter
- inability to push the retrieval catheter over the filter
- restriction of blood flow
- distal embolization
- infection
- filter fracture.

Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.

Retrieval of the OPTEASE Retrievable Filter with a filter fracture present may result in complications.

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1. Prior to commencing the retrieval procedure, determine that the retrieval route is free of thrombus.
2. Verify the duration of the OPTEASE Retrieval Filter implantation. Do not retrieve the OPTEASE Retrieval Filter if it has been implanted in the patient for longer than 12 days.
3. Remove the OPTEASE Retrieval Catheter from its packaging using sterile technique.
4. Inspect the OPTEASE Retrieval Catheter upon removal from packaging to verify that it is undamaged. Warning: Do not use an OPTEASE Retrieval Catheter that has been damaged in any way. If damage is detected, replace with an undamaged OPTEASE Retrieval Catheter.
5. Prior to use, flush the OPTEASE Retrieval Catheter with heparinized saline or suitable isotonic solution.
6. Prepare all other procedure components according to the manufacturers' *Instructions for Use*.
7. Attach the hemostasis device to the hub of the OPTEASE Retrieval Catheter.
8. Use appropriate techniques to determine that the filter, the femoral retrieval route, and distal IVC are free of thrombus.
9. Select the correct loop diameter size of the Microvena Amplatz Goose Neck Snare as indicated in Table 1.
10. Assemble the Microvena Amplatz Goose Neck Snare Kit according to the *Instructions for Use*.
11. Using the guidewire tip straightener, insert the recommended guidewire into the 10F catheter sheath introducer. Gently advance the guidewire in the IVC such that it is caudal to, but not into, the OPTEASE Retrieval Filter.

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12. Introduce and advance the OPTEASE Retrieval Catheter such that the tip of the OPTEASE Retrieval Catheter is located at a short distance (approximately 3 cm) caudal of the OPTEASE Retrieval Filter retrieval hook.

13. Remove the guidewire.

14. Insert and advance the Microvena Amplatz Goose Neck Snare assembly through the OPTEASE Retrieval Catheter until it protrudes out of the OPTEASE Retrieval Catheter such that the marker band of the snare catheter is caudal of the filter retrieval hook.

15. Push the snare shaft forward gently to open the snare loop caudal of the OPTEASE Retrieval Filter retrieval hook.

16. The loop is then slowly advanced forward over the OPTEASE Retrieval Filter apex (Figure 2).

17. Reduce the loop diameter by advancing the snare catheter while simultaneously pulling backwards the snare until the loop engages the filter retrieval hook (Figure 3A). Note: Ensure that the loop of the snare has properly engaged the retrieval hook and that the retrieval hook, retrieval catheter and snare are aligned (Figure 3A illustrates correctly aligned, Figure 3B illustrates incorrectly aligned). The marker tip of the snare catheter must be caudal to the filter hook. There must be a clear distance between the marker tip of the Microvena Snare Catheter and the denser end of the OPTEASE Retrieval Filter apex.

18. Always maintain tension on the snare to prevent disengagement of the snare loop from the OPTEASE Retrieval Filter retrieval hook.

19. While keeping tension on the snare, ensure that the filter and snare catheter are in line with the OPTEASE Retrieval Catheter (Figure 3A). Pull the snare catheter approximately 5 mm backward to allow the loop to align with the filter.

20. While keeping tension of the snare, advance the OPTEASE Retrieval Catheter over the filter hook (Figure 4). Continue advancing the OPTEASE Retrieval Catheter over the filter until the filter is completely collapsed inside the OPTEASE Retrieval Catheter. Precaution: Do not push the snare catheter against the retrieval hook prior to the advance of the OPTEASE Retrieval Catheter over the filter. Too close contact between the snare catheter and

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filter retrieval hook can cause friction of the filter tip in the OPTEASE Retrieval Catheter. Note: If for any reason the OPTEASE Retrieval Filter is not retrieved and remains implanted as a permanent filter, perform fluoroscopic control to determine the status of the IVC and filter before terminating the procedure.

21. Once the filter is fully collapsed inside the OPTEASE Retrieval Catheter retract the complete system as a unit out through the 10F sheath introducer. Warning: DO NOT retract the filter out of the OPTEASE Retrieval Catheter through the valve of the hemostasis device. This could damage the valve of the hemostasis device.

22. Control the status of the IVC before ending the procedure using an appropriate technique.

23. Remove the sheath introducer when clinically indicated by placing compression on the vessel above the puncture site, and slowly withdrawing the system.

24. Discard the OPTEASE Retrieval Filter, OPTEASE Retrieval Catheter and packaging materials. Note: After use, all components used and packaging materials may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable local, state and federal laws and regulations.

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able Filter should only be performed by physicians who are trained in percutaneous interventional techniques. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.

The use of any other components could result in complications and/or an unsuccessful procedure. Cordis recommends the use of the Microvena Amplatz Goose Neck Snare Kit for the retrieval of the OPTEASE Retrievable Filter. The OPTEASE Retrievable Filter and OPTEASE Retrieval Catheter have been tested with the Microvena Amplatz Goose Neck Snare Kit.

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The guiding OUTBACK LTD Re-Entry Catheter (OB-LTD) is designed to facilitate placement and positioning of guidewires within the peripheral vasculature. The OB-LTD is comprised of three primary elements: 1) Cannula, 2) Catheter shaft, and 3) Deployment handle with deployment control slide.

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The OUTBACK LTD Re-Entry guiding Catheter is intended to facilitate placement and positioning of guidewires and catheters within the peripheral vasculature.

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The OUTBACK LTD Re-Entry Catheter is not intended for use in the coronary or cerebral vasculature.

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This product is designed for use by physicians trained in and familiar with peripheral percutaneous interventional techniques.

Possible complications may include, but are not limited to, the following:

- Vessel dissection, perforation, or injury
- Embolism
- Vessel spasm
- Pseudoaneurysm
- Vascular thrombosis

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

This medical device is used only in the conditions of medical and therapeutic in-

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stitutions by qualified medical personnel.

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the handle deployment slide is locked in the most proximal position. Always track the OB-LTD over a guide wire. Do not track the OB-LTD in the vasculature without a guide wire.

9. Track the OB-LTD over the guide wire to the desired vascular site. Torque the OB-LTD as needed during delivery via the handle RHV.
10. Retract the guide wire tip into the catheter approximately 5 cm, using fluoroscopic guidance to confirm guide wire position. Failure to retract the guide wire tip into the catheter approximately 5 cm before deployment of cannula could result in fracture or separation of wire tip from body of wire.
11. Rotate the image intensifier so the distal housing of the OB-LTD is adjacent to the target re-entry site when visualized using fluoroscopy.
12. Orientate the lateral exit port of the OB-LTD towards the desired vascular target site via rotation of the RHV. Using fluoroscopic guidance, position the "L" marker leg, on the catheter "LT" directional marker band, towards the target re-entry site using an initial fluoroscopic view (See Fig 1)
13. Once this initial orientation is confirmed, rotate the image intensifier to a ninety-degree orthogonal fluoroscopic view and confirm that the distal housing of the OB-LTD is positioned "in line" with target re-entry site.
14. Orientate the lateral exit port of the OB-LTD towards the desired vascular target site via rotation of the RHV. Using fluoroscopic guidance, rotate RHV until the catheter "LT" directional marker band on the nosecone appears as a "T". See Fig 2

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ing the catheter over the guide wire

21. Carefully retract the OB-LTD over the guide wire, leaving the guide wire in place for subsequent therapeutic procedures

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The following guide wires are recommended for use with the OUTBACK LTD Re-Entry Catheter (OB-LTD). LuMend test results have shown these wires to be compatible with the OB-LTD. Failure to use a recommended guide wire may result in damage to the guide wire, such as, abrasion of the hydrophilic coating, release of polymer fragments, separation of the wire, or inability to withdraw the OB-LTD over the guide wire.

.014" ATW (Cordis, Johnson & Johnson Company)
.014" STABILIZER® Plus (Cordis, Johnson & Johnson Company)
.014" STABILIZER® XS (Cordis, Johnson & Johnson Company)
.014" Sparta Core (Guidant Corporation)
.014" Iron Man (Guidant Corporation)
.014" Whisper (Guidant Corporation)
.014" All Star (Guidant Corporation)
.014" Platinum Plus (Boston Scientific/Scimed)
.014" Choice Extra Support (Boston Scientific/Scimed)
.014" Mailman (Boston Scientific/Scimed)
.014" Luge (Boston Scientific/Scimed)

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aged.

- Do not use after the last day of the month of the "Use By" date on the package

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proximal hub to the distal tip. Do not tug or otherwise overstretch the catheter to straighten it.

- Always confirm visualization of the targeted distal vessel via contrast injection and fluoroscopy • before using the catheter. Avoid contrast injection in the sub-intimal space.
- Minimize sub-intimal dissection tract beyond point of reconstitution.
- To maintain guide wire position during device exchanges, an exchangeable length guide wire is recommended.
- Prior to use, carefully read the instructions packaged with the guide wire that is to be used with the OB-LTD. See Table 1 for a list of recommended guide wires. Failure to use a recommended guide wire may result in damage to the guide wire, such as abrasion of the hydrophilic coating, release of polymer fragments, separation of the wire, or inability to withdraw the OB-LTD over the guide wire.
- Always visualize tracking of the catheter tip over the aorto-iliac bifurcation.
- If strong resistance is felt during catheter manipulation/delivery, determine the cause of the resistance before proceeding further. Consider using a 3–4 mm balloon at low ATM to dilate points of resistance, as needed, along delivery track. If the cause cannot be determined, withdraw the OB-LTD.
- Excessive rotation, bending or kinking of the OB-LTD may affect its performance. Withdraw the OB-LTD if it becomes excessively kinked.

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T (786) 313 2000 F (786) 313 2080

April 20, 2017

CERTIFICATION OF TRUTHFUL AND ACCURATE STATEMENT

I, **Miguel Ávila**, in my capacity as Vice President, Global Quality and Regulatory Affairs – Cordis Corporation certify that the attached **Instructions for Use for Guidewires**:

- ATW;
- ATW Marker wire;
- EMERALD;
- JINDO;
- REFLEX;
- SHINOBI;
- SHINOBI PLUS;
- STABILIZER MARKER WIRE;
- STABILIZER PLUS;
- STABILIZER EXTRA SUPPORT;
- STABILIZER SUPPORT;
- STABILIZER BALANCED PERFORMANCE;
- STABILIZER XS;
- STORQ;
- SV;
- WIZDOM;
- CINCH QR;
- ANGIOGUARD RX;
- FrontRunner XP;

issued by **Cordis Corporation** are true and correct to the best of my knowledge.

Miguel Ávila

Vice President, Global Quality and Regulatory Affairs


State of California
County of Santa Clara

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document, to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

On April 20, 2017, before me, ZHENG WANG, Notary Public, personally appeared Miguel Avila, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity, and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

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- ATW
- SHINOBI
- SHINOBI PLUS
- STABILIZER MARKER WIRE
- STABILIZER PLUS
- STABILIZER EXTRA SUPPORT
- STABILIZER SUPPORT
- STABILIZER BALANCED PERFORMANCE
- STABILIZER XS
- WIZDOM

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Cordis guidewires are intended for use in angiographic procedures to introduce and position catheters and interventional devices within the coronary and peripheral vasculature.

Cordis steerable guidewires are contraindicated for use in chronic total occlusions.

Contraindications for interventional devices are described in the instructions supplied with the respective device.

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- air embolism
- hematoma at the puncture site
- infection

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- perforation of the heart or vessel wall
- Tip fractures have been reported in procedures involving guidewire entrapment, total occlusions, highly tortuous vasculature, and small side branches

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1. Use appropriately sized Steerable Guidewires with the selected interventional device (refer to the interventional device for compatible guidewire size).

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3. The guidewire can be loaded into the interventional device by using one of two methods (refer to the interventional device's instructions for use for the preferred

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method):

а) Insert the proximal end of the steerable guidewire into the interventional device tip and advance the guidewire until it exits the guidewire hub.

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Caution: To prevent vessel injury or **tip entrapment**, use fluoroscopy when advancing/torquing the guidewire.

5. Use accepted angiographic technique for guidewire positioning. Then continue the procedure.

Caution: If strong resistance is met during manipulation, discontinue the procedure and determine cause of resistance before proceeding. If the cause of resistance can not be determined, withdraw the entire system. If the tip becomes fixed within the vasculature, advance the balloon catheter as distally as possible, gently pull the guidewire back into the balloon catheter and withdraw the balloon catheter/guidewire system as a unit - **NEVER TORQUE THE GUIDEWIRE IF THE TIP BECOMES FIXED.**

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Note: When the guidewire is in the interventional device's inner lumen, pressure

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measurement and dye injection through the inner lumen may be impaired. Consult the device's instructions for the recommended procedure in these circumstances.

7. To change devices, withdraw the interventional device into the guiding catheter while maintaining guidewire position within the vasculature.

Continue to withdraw the device until the distal tip exits the guide catheter's proximal hub.

Remove the device and backload the second device over the guidewire and repeat steps 2 through 5.

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- Guidewire manipulation/torquing should

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always be performed under fluoroscopic guidance.

- Never push, auger, withdraw, or torque a guidewire that meets resistance. First, using fluoroscopy, determine the cause of resistance and take any necessary remedial action. Torquing or pushing a guidewire against resistance may cause guidewire damage, and/or guidewire tip separation, or direct damage to the vessel. Resistance may be felt and/or observed (via fluoroscopy) by noting any buckling of the guidewire tip. If guidewire tip prolapse is observed, **DO NOT** allow the tip to remain in a prolapsed position; otherwise damage to the guidewire may occur.

- If any resistance is felt, i.e., due to vessel spasm, bent guidewire, or guidewire entrapment, while manipulating or removing the guidewire in the blood vessel: **STOP** the procedure. **DO NOT** move or torque the guidewire. Using fluoroscopy, first determine the cause of the resistance, then take appropriate remedial action. If the guidewire is moved excessively, it may break or become damaged. This may cause blood vessel injury or result in fragments being left inside the vessel.

- Should torque control/tip response be compromised during use, confirm tip integrity using fluoroscopy. **LOSS OF TORQUE CONTROL MAY BE DUE TO CORE WIRE FRACTURE.** Under fluoroscopic guidance, advance the balloon catheter to the distal end of the guidewire and remove the balloon catheter/guidewire system as a unit.

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Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the

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product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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Instruction for Use Guidewires

The Cordis Steerable Guidewires have a soft, flexible, highly radiopaque platinum alloy coil on the distal tip. Guidewire length, diameter, and distal tip configuration are indicated on the product label. A steering/torquing device may be attached to aid in steering/torquing the guidewire.

Cordis guidewires are intended for use in angiographic procedures to introduce and position catheters and interventional devices within the coronary and peripheral vasculature. In addition, the guidewires are intended to facilitate the alignment of interventional devices and function as a measurement tool.

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Procedures requiring percutaneous guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at anytime during or after the procedure.

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- air embolism
- hematoma at the puncture site
- infection
- perforation of the heart or vessel wall
- Tip fractures have been reported in procedures involving guidewire entrapment, total occlusions, highly tortuous vasculature, and small side branches..

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1. Use appropriately sized steerable guidewires with the selected interventional device (refer to the interventional device for compatible guidewire size).

2. Open the sterile package slowly. To prevent damage to the distal tip of the guidewire and to prevent the guidewire from springing onto a non-sterile field, carefully remove the guidewire from the dispensing tube. Use the dispenser window to advance the guidewire distal tip, such that the guidewire can be removed by grasping the coated shaft.

Note: Do not pull the distal tip to remove guidewire from dispenser as it may damage the tip.

Caution: Before use, flush all devices entering the vascular system with sterile, heparinized saline or a similar isotonic solution.

3. The guidewire can be loaded into the interventional device by using one of two methods (refer to the interventional device's Instructions for Use for the preferred method):

a) insert the proximal end of the steerable guidewire into the interventional device tip and advance the guidewire until it exits the

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guidewire hub of the interventional device, or

b) insert the distal tip of the guidewire into the guidewire hub of the interventional device and advance the guidewire until the tip exits the interventional device's tip.

Caution: Gently introduce and advance the guidewire to prevent damaging the distal tip.

Note: To aid in the navigation within and selection of arteries, the guidewire may be reshaped.

4. To aid in rotating or steering the guidewire, secure a steering device to the proximal end of the guidewire.

Caution: To prevent vessel injury or tip entrapment, use fluoroscopy when advancing/torquing the guidewire.

5. Use accepted angiographic technique for guidewire positioning.

Then continue the procedure.

Caution: If strong resistance is met during manipulation, discontinue the procedure and determine cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system. If the tip becomes fixed within the vasculature, advance the balloon catheter as distally as possible, gently pull the guidewire back into the balloon catheter and withdraw the balloon catheter/ guidewire system as a unit - **NEVER TORQUE THE GUIDEWIRE IF THE TIP BECOMES FIXED.**

6. When desired procedure results are achieved, withdraw the guidewire and interventional device slowly.

Note: When the guidewire is in the interventional device's inner lumen, pressure measurement and dye injection through the inner lumen may be impaired. Consult the device's instructions for the recommended procedure in these circumstances.

7. To change devices, withdraw the interventional device into the guiding catheter while maintaining guidewire position within the vasculature. Continue to withdraw the device until the distal tip exits the guide catheter's proximal hub. Remove the device and backload the second device over the guidewire and repeat steps 2 through 5.

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Caution: Device removal should be performed under fluoroscopic guidance.

Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.

Guidewires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends, or kinks. Do not use a guidewire that shows signs of damage. Damage will prevent the guidewire from performing with accurate torque response and control.

- Guidewire manipulation/torquing should always be performed under fluoroscopic guidance.
- Never push, auger, withdraw, or torque a guidewire that meets resistance. First, using fluoroscopy, determine the cause of re-

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сistance and take any necessary remedial action. Torquing or pushing a guidewire against resistance may cause guidewire damage, and/or guidewire tip separation, or direct damage to the vessel. Resistance may be felt and/or observed

(via fluoroscopy) by noting any buckling of the guidewire tip.

If guidewire tip prolapse is observed, DO NOT allow the tip

to remain in a prolapsed position; otherwise damage to the guidewire may occur.

- If any resistance is felt, e.g., due to vessel spasm, bent guidewire, or guidewire entrapment, while manipulating or removing the guidewire in the blood vessel: STOP the procedure. DO NOT move or torque the guidewire. Using fluoroscopy, first determine the cause of the resistance, then take appropriate remedial action. If the guidewire is moved excessively, it may break or become damaged. This may cause blood vessel injury or result in fragments being left inside the vessel.

- Should torque control/tip response be compromised during use, confirm tip integrity using fluoroscopy. LOSS OF TORQUE CONTROL MAY BE DUE TO CORE WIRE FRACTURE. Under fluoroscopic guidance, advance the balloon catheter to the distal end of the guidewire and remove the balloon catheter/guidewire system as a unit.

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Descriptions or specifications in Cordis

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printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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Sterility maintaining term - 2 years.

After use all used components and packing materials may be potentially biologically dangerous. Their treatment and disposal are performed according to the accepted medical practice, and also according to local, state and federal laws and provisions.

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Instruction for Use Guidewires EMERALD

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Cordis Guidewires have a stainless steel shaft. The shaft is either PTFE* coated, PTFE and Heparin coated, or uncoated. Guidewire length, diameter, core configuration, flexibility type, coating, and tip configuration are all indicated on the product label.

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Cordis Guidewires are intended for use in the percutaneous introduction of catheters.

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Procedures requiring percutaneous catheter/guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

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1. Flush the guidewire with sterile heparinized normal saline or a similar isotonic solution by connecting a syringe to the Luer hub of the guidewire dispenser.
2. Insert the flexible end of the guidewire into the needle.
3. Advance the guidewire through the needle.
4. Position the guidewire.
5. Hold the guidewire in place and withdraw the needle.
6. Advance the catheter over the guidewire.

Movable Core Guidewires

To lengthen the flexible segment of a movable core guidewire, gently withdraw the core. To shorten the flexible segment, push the core in carefully.

Caution: When the guidewire is in a vessel, do not advance the movable core if the tip is in a curved shape. Never twist or force the core because excessive force may cause it to penetrate the coil and damage the vessel.

Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Accordingly, Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse

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Do not reuse. Discard after one procedure. Structural integrity and/ or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.

Guidewires are delicate instruments and should be handled carefully.

Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends, or kinks. Do not use a guidewire that shows signs of damage. Damage will prevent the guidewire from performing with accurate torque response and control.

- Guidewire manipulation/torquing should always be performed under fluoroscopic guidance.

- Never advance, withdraw or auger the guidewire against resistance without first determining the cause of resistance under

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fluoroscopy. Torquing the guidewire against resistance may cause damage and/or fracture which may result in separation of the distal tip.

- Should the guidewire tip become entrapped within the vasculature (i.e., small side branch, tight stenosis), **DO NOT TORQUE THE GUIDEWIRE.** Advance the balloon catheter distally, gently pull the guidewire back into the balloon catheter, and remove the balloon catheter/guidewire system as a unit.

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- Store in a cool, dark, dry place.
- Do not use open or damaged packages.

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After use all used components and packing materials can be potentially biologically dangerous. Disposal is carried out according to the accepted medical practice, and also according to local, state and federal laws and provisions.

The products have come into contact with blood and / or body fluids, disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class).

Unused products (which had no contact

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with blood and / or body fluids), including expired, disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

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Instruction for Use Guidewires REFLEX

The Cordis Steerable Guidewires have a soft, flexible, highly radiopaque platinum/tungsten coil on the distal tip. Guidewire length, diameter, and distal tip configuration are indicated on the product label. A steering/torquing device may be attached to aid in steering/torquing the guidewire.

Cordis Steerable Guidewires are intended for use in angiographic procedures to introduce and position catheters and interventional devices within the vasculature.

Cordis Steerable Guidewires are contraindicated for use in chronic total occlusions. Contraindications for interventional devices are described in the instructions supplied with the respective device.

Procedures requiring percutaneous guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at anytime during or after the procedure. Possible complications include, but are not limited to:

- Air embolism
- Hematoma at the puncture site
- Infection
- Perforation of the heart or vessel wall
- Tip fractures have been reported in proce-

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dures involving guidewire entrapment, total occlusions, highly tortuous vasculature, and small side branches.1-6 For guidewire tip retrieval, please refer to the referenced publications for recommended techniques.

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3. The guidewire can be loaded into the interventional device by using one of two methods (refer to the interventional device's instructions for use for the preferred method):

a) Insert the proximal end of the steerable guidewire into the interventional device tip and advance the guidewire until it exits the guidewire hub.

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б) Insert the distal tip of the guidewire into the guidewire hub of the interventional device and advance the guidewire until the tip exits the interventional device's tip.

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Then continue the procedure.

Caution: If strong resistance is met during manipulation, discontinue the procedure and determine cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system. If the tip becomes fixed within the vasculature, advance the balloon catheter as distally as possible, gently pull the guidewire back into the balloon catheter and withdraw the balloon catheter/guidewire system as a unit - **NEVER TORQUE THE GUIDEWIRE IF THE TIP BECOMES FIXED.**

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- Never push, auger, withdraw, or torque a guidewire that meets resistance. First, using fluoroscopy, determine the cause of resistance and take any necessary remedial action.

Torquing or pushing a guidewire against resistance may cause guidewire damage, and/or guidewire tip separation, or direct damage to the vessel. Resistance may be felt and/or observed (via fluoroscopy) by noting any buckling of the guidewire tip. If guidewire tip prolapse is observed, DO NOT allow the tip to remain in a prolapsed position; otherwise damage to the guidewire may occur.

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Instruction for Use Guidewires

The Cordis Steerable Guidewire Extension is a flexible stainless steel wire coated with DURAGLIDE coating. The extension is approximately 145 cm long and has a shaft diameter of 0.014 inches

The Cordis Steerable Guidewire Extension is designed to extend Cordis Steerable Guidewires so that a catheter can be exchanged for another catheter.

None known.

Procedures requiring percutaneous guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure. Possible complications include, but are not limited to the following:

- air embolism
- hematoma at puncture site
- infection
- perforation of the heart or vessel wall

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1. Open the outer box and remove the pouch. Open the sterile pouch slowly and remove the dispensing tube containing the extension.

2. Compress the winged alignment tool between the thumb and forefinger. This holds the guidewire extension stationary.

3. Place the proximal end of the Cordis guidewire into the funneled end of the blue winged alignment tool.

4. Gently push the guidewire until fully seated within the extension wire.

Note: Do not attempt to use or straighten an extension which exhibits damage to the hypotube. Attempts to straighten the hypotube may result in fracturing of the hypotube.

Verify that proper engagement has occurred by gently pulling on the guidewire. If it appears to spring back towards the extension proceed to Step 5. If it appears to be sliding out of the extension, then verify that the guidewire is actually in the extension hypotube and gently reinsert the guidewire into the extension. Repeat as required.

5. Carefully lift and pull the extended guidewire out of the winged alignment tool and remove from the dispensing tube.

Caution: Before use, flush all devices entering the vascular system with sterile, heparinized saline or similar isotonic solution.

Note: To minimize the potential of inadvertent disengagement during wipe down, it is recommended that the guidewire/extension connection be firmly grasped during the wipe down.

6. Exchange using standard procedures.

7. Disengage the extension from the guidewire by one of the following methods:

a. The recommended procedure is to locate the extension's hypotube visually, and/or by tactile means by running fingers along the steerable guidewire in the proximal direction. While holding the guidewire firmly in place, move the extension device's hypotube away from the guidewire until it stops. While holding the hypotube in that posi-

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тион, disconnect the extension device by pulling the extension away from the guidewire.

б. An alternative method for disengaging the guidewire from the extension is to locate the extension's hypotube. Grip the extension's wire just proximal to the hypotube and grip the guidewire just distal to the hypotube. Rotate the extension device counterclockwise and simultaneously pull it away from the steerable guidewire until they disconnect.

Do not reuse. Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Accordingly, Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.

Never push, auger, withdraw, or torque a guidewire that meets resistance. **STOP** the procedure. **DO NOT** move or torque the guidewire. First, using fluoroscopy, determine the cause of resistance and take any necessary

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remedial action. Torquing or pushing a guidewire against resistance may cause guidewire damage, and/or guidewire tip separation, or direct damage to the vessel. Resistance may be felt and/or observed (via fluoroscopy) by noting any buckling of the guidewire tip. If guidewire tip prolapse is observed, DO NOT allow the tip to remain in a prolapsed position; otherwise damage to the guidewire may occur.

Guidewires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guidewire carefully for signs of damage. Do not use a guidewire that shows signs of damage.

- Store in a cool, dark, dry place.
- Do not use open or damaged packages.
- Use prior to "Use By" date.
- Do not resterilize. Exposure to temperatures above 54°C (130°F) may damage the extension.
- Do not expose to organic solvents.
- If the hypotube at the distal end of the guidewire extension is damaged, do not use. Do not attempt to straighten a kinked hypotube. Replace the unit with an undamaged unit before proceeding.

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Instruction for Use Guidewires ANGIOGUARD RX

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The ANGIOGUARD RX Emboli Capture Guidewire is a steerable guidewire available in several filter basket diameters. Guidewire length, diameter, distal tip configuration and guidewire support are indicated on the product label. It is a single use, disposable guidewire with the addition of a filter basket at the distal end of the guidewire. The filter basket is used for trapping emboli during coronary, carotid, and peripheral procedures. It consists of a thin, porous membrane supported by a fine metal skeleton. In the collapsed state, the filter basket has a very low profile. The ANGIOGUARD RX Emboli Capture Guidewire System is used in a standard fashion to cross the lesion and support the placement of devices such as angioplasty balloons and intravascular stents. Once the guidewire is across the lesion, the filter basket is expanded in an umbrella-like fashion in the vessel lumen. During the coronary, carotid or peripheral procedure, as the blood passes through the filter basket, emboli are collected. At the completion of the procedure, the filter basket is collapsed, capturing the emboli, so it can be removed through the guiding catheter or interventional sheath Introducer and out of the body. The ANGIOGUARD RX Emboli Capture Guidewire System has two (2) exit markers located on the proximal shaft of the deployment and capture sheaths. These markers, located at 90 and 100 cm respectively, approximate the exit of the guidewire distal tip from the guiding catheter or the interventional sheath.

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The Cordis ANGIOGUARD RX Emboli Capture Guidewire System consists of one (1) each: ANGIOGUARD RX Emboli Capture Guidewire, deployment sheath, capture sheath, torque device/locking device, filter basket introducer, and peel away guidewire introducer. (See Figure 1.)

The ANGIOGUARD RX Emboli Capture Guidewire System is intended for coronary, carotid, and peripheral uses to facilitate the placement of diagnostic and interventional devices, and capture emboli, thereby reducing the risk of embolization, during procedures.

The ANGIOGUARD RX Emboli Capture Guidewire System is contraindicated for use in conjunction with devices to treat intracranial lesions. The ANGIOGUARD RX Emboli Capture Guidewire System is contraindicated for use in chronic total occlusions.

Procedures requiring percutaneous guidewire introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

Possible complications include, but are not limited to, the following:

Coronary Indication:

- Air embolism
- Infection
- Vessel dissection, perforation, rupture, or injury
- Death

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- Acute myocardial infarction
- Unstable angina
- Dysrhythmia, including ventricular fibrillation
- Hypo/Hypertension
- Vessel Spasm
- Arteriovenous or vessel fistula
- Intravascular stent migration

Carotid Indication:

- Air embolism
- Infection
- Vessel dissection, perforation, rupture, or injury
- Death
- Acute myocardial infarction
- Unstable angina
- Dysrhythmia, including ventricular fibrillation
- Hypo/Hypertension
- Vessel Spasm
- Arteriovenous or vessel fistula
- Intravascular stent migration
- TIA
- Stroke

Peripheral indication:

- Air embolism
- Infection
- Vessel dissection, perforation, rupture, or injury
- Vessel spasm
- Arteriovenous or vessel fistula
- Ischemia
- Intravascular stent migration
- TIA
- Stroke

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Do not use the device if there are abnormalities in the product or in the sterile barrier (e.g. broken seal, torn or breached barrier).

This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

Use the ANGIOGUARD RX Emboli Capture Guidewire System prior to the "Use By" date specified on the package.

Observe all ANGIOGUARD RX Emboli Capture Guidewire System movement in the vessels using fluoroscopic guidance. Before the guidewire is moved, the tip should be visible using fluoroscopy. Do not torque a guidewire without observing corresponding movement of the tip; otherwise vessel trauma could occur.

If the wire tip becomes entrapped within the vasculature, DO NOT TORQUE THE GUIDEWIRE.

Torquing a guidewire against resistance may cause guidewire damage and/or guidewire tip separation. Always advance or withdraw the guidewire slowly. Never push, auger, withdraw or torque a guidewire, which meets resistance. Resistance may be felt and/or observed using fluoroscopy by noting any buckling of the guidewire tip. Determine the cause of resistance using fluoroscopy and take the necessary remedial action.

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Perform all exchanges slowly to prevent air from entering the catheter system or trauma.

When introducing the guidewire, confirm that the guiding catheter or interventional sheath Introducer tip is free within the vessel lumen and not against the vessel wall. Failure to do so may result in vessel trauma upon guidewire exit from the device tip. Use the radiopaque marker of the interventional device (guiding catheter or sheath) to confirm position.

Maintain an ACT of > 300 seconds during the time the ANGIOGUARD RX Emboli Capture Guidewire filter is deployed.

Animal studies demonstrated the patency of this device *in vivo* for up to 60 minutes.

- If distal perfusion of dye is significantly reduced or no dye is perfusing past the distal marker band, the ANGIOGUARD RX Emboli Capture Guidewire filter may have reached its maximum capacity to contain emboli. Remove and replace with a new ANGIOGUARD RX Emboli Capture Guidewire System per corresponding section.

- Do not attempt to close the filter basket with the deployment sheath. The ANGIOGUARD RX Emboli Capture Guidewire should only be removed using the capture sheath.

- Care during diagnostic or interventional device exchanges must be practiced to minimize movement of the guidewire filter basket.

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- Guidewires are delicate instruments and should be handled carefully. Prior to use and when possible during the procedure, inspect the guidewire carefully for coil separation, bends, kinks, or damage of the filter basket assembly. Figures 1 and 2 are used for reference only to show coil dispenser, torque device, filter basket introducer, emboli capture guidewire and basket, deployment sheath, capture sheath, and flushing Luer hub.

Figure 1



Figure 2

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Preparations for Use

Prior to the interventional procedure, all equipment and packaging, including the ANGIOGUARD RX Emboli Capture Guidewire System, should be inspected and examined carefully for defects. Check guidewire, filter basket, deployment sheath, capture sheath, and capture sheath RX port region (approximately 30 cm from the distal tip) for bends, kinks, or other damage. Do not use defective equipment.

1. The package contains two plastic coil dispensers. One contains the guidewire with filter basket, filter basket introducer, peel away guidewire Introducer, torque device, and black/yellow deployment sheath. The second contains the blue capture sheath, which will not be needed until the filter basket is ready to be removed from the vessel.

2. During shipping the deployment sheath tip may become disengaged from the filter basket introducer. Verify deployment sheath tip is engaged. If not, engage manually by inserting deployment sheath tip into filter basket introducer.

3. Fill a 10-ml Luer lock syringe with sterile saline and purge all air from syringe.

4. Attach syringe to Luer lock hub on the end of the filter basket introducer.

5. Inject 10 ml of saline to purge all air from the deployment sheath and filter basket (ensure distal tip of the deployment sheath is engaged with the filter basket introducer tip prior to purging the system). You may see saline dripping along the length of the deployment sheath. Flushing is complete when saline is seen within coil dispenser at the green deployment sheath hub.

6. Disconnect syringe.

7. Remove the two anti-migration clips closest to torque device and pull to ensure torque device is securely attached to guidewire.

8. Gripping the torque device in one hand and the coil dispenser in the other, pull on the wire until the basket is completely docked into the tip of the deployment sheath. When completely docked, approximately half the filter basket will still be visible out the end

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of the deployment sheath.

9. Remove the last anti-migration clip.

10. Open the torque device.

11. Gripping the torque device in one hand and the proximal end of the guidewire in the other, pull the wire through the torque device until the proximal end of the deployment sheath hub engages the torque nut.

12. Lock the torque device onto the guidewire; ensure the deployment sheath hub remains engaged with the torque nut. Pull the wire and deployment sheath out of dispenser coil.

13. The deployment sheath is now prepped and ready for use.

14. Capture sheath – After the interventional or diagnostic procedure is completed, attach a 10-ml Luer lock syringe to flushing Luer hub, prepared as described in Preparations for Use, step 3. Flush with 10 ml of sterile saline. The entire coil dispenser should be filled with saline. Remove the syringe from the flushing Luer hub. Do not remove the capture sheath from the coil dispenser until required system retrieval.

Filter Basket Sizing Recommendation

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Nominal Filter Basket Size
4.0 mm
5.0 mm
6.0 mm
7.0 mm
8.0 mm

Recommended Vessel Size
3.0 mm < Vessel ≤ 3.5 mm
3.5 mm < Vessel ≤ 4.5 mm
4.5 mm < Vessel ≤ 5.5 mm
5.5 mm < Vessel ≤ 6.5 mm
6.5 mm < Vessel ≤ 7.5 mm

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Directions for use

1. Insert the peel away guidewire introducer into the interventional sheath introducer or the hemostatic valve of the Y connector attached to the guiding catheter. Ensure the hemostatic valve is fully opened before inserting the guidewire intro-

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

nique.

Figure 3

11. The ANGIOGUARD RX Emboli Capture Guidewire is now fully deployed within the vessel. Confirm full deployment using fluoroscopy. The marker bands on the filter basket struts should be fully opposed to the vessel wall.
12. Using fluoroscopy, inject dye to determine that there is adequate flow distal to the filter basket (or distal band of the guidewire) and that the guidewire remains in the proper position. Note: Allow sufficient space between the lesion and filter basket to prevent interference with the distal tips of all other diagnostic and interventional devices (i.e. stent delivery systems, angioplasty balloons) to be used throughout the procedure.
13. Once guidewire position is confirmed, any rapid exchange .014" (0.36 mm) compatible interventional device can be loaded onto the guidewire for treatment of the lesion. Use the standard rapid exchange technique while treating the lesion. Be careful not to move the filter basket during exchanges. Note: Never attempt to capture the filter basket with the deployment sheath.
14. Once the lesion is treated and all of the interventional or diagnostic devices have been removed, remove the flushed capture sheath (per step

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14 of section IX) from the dispenser coil and thread over the proximal end of the guidewire. Grasp the guidewire proximally as it exits the RX port and push the capture sheath through the open hemostasis valve. Collapse the filter basket by advancing the capture sheath until the distal capture sheath marker band is adjacent to the proximal filter basket marker band. Using fluoroscopy, confirm filter basket closure by ensuring reduction of diameter of the radiopaque strut marker bands. Note: Do not try removing the ANGIOGUARD RX Emboli Capture Guidewire by only pulling on the capture sheath. Refer to Figure 4. Note: Never pull the ANGIOGUARD RX Emboli Capture Guidewire into the guiding catheter or interventional sheath Introducer if there is resistance. If resistance is encountered, reposition the capture sheath to ensure that the filter basket is properly seated into the capture sheath. Using fluoroscopy, verify capture sheath position by checking marker band alignment between the capture sheath and the guidewire.

15. Remove the device by grasping the guidewire and capture sheath together (refer to Figure 4) near the hemovalve. Pull the system through the guiding catheter or interventional sheath Introducer, and out of the hemovalve as a single unit. Care should be taken when pulling the basket through the open hemovalve, to avoid potential release of captured emboli.

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

Figure 4

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How to Detect if the ANGIOGUARD® RX Emboli Capture Guidewire Is Full

1. Once the ANGIOGUARD RX Emboli Capture Guidewire is deployed in the vessel, it may capture emboli during the entire time of the interventional procedure. Therefore, it is recommended to check the status of the ANGIOGUARD RX Emboli Capture Guidewire at regular intervals during the intervention.

2. Using fluoroscopy, perform a distal dye injection through the guiding catheter or interventional sheath Introducer and observe the flow of dye distal to the filter basket or distal marker on the guidewire.

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3. If the distal perfusion of dye is significantly reduced or no dye is perfusing past the filter basket or guidewire distal marker band, the ANGIOGUARD RX Emboli Capture Guidewire may have reached its capacity to contain emboli. If there is a severe reduction in distal dye perfusion, it is recommended to exchange the ANGIOGUARD RX Emboli Capture Guidewire for a new one.

How to Exchange a Full ANGIOGUARD RX Emboli Capture Guidewire

Caution: The ANGIOGUARD RX Emboli Capture Guidewire should only be removed using the capture sheath and should be removed through the guiding catheter or interventional sheath Introducer ONLY. See Section X, steps 14 & 15 for proper filter basket capture and system removal.

To exchange an ANGIOGUARD RX Emboli Capture Guidewire System that has captured its capacity of emboli:

1. Remove all interventional devices from the ANGIOGUARD RX Emboli Capture Guidewire.
2. Prep the capture sheath as outlined in the Preparations for Use section, step 14 of section IX.
3. Load the capture sheath over the proximal end of the ANGIOGUARD RX Emboli Capture Guidewire. Grasp the guidewire proximally as it exits the RX port and advance the capture sheath through the open hemostasis valve. Continue to advance the capture sheath until the distal capture sheath marker band is adjacent to the proximal filter basket marker band. This collapses the filter basket. Using fluoroscopy, confirm filter basket closure by ensuring reduction of diameter of the radiopaque strut marker bands. **Note:** The emboli that have been captured may not allow the ANGIOGUARD RX Emboli Capture Guidewire to reach its initial low profile.

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4. Open the hemostasis valve to allow the ANGIOGUARD RX Emboli Capture Guidewire free movement and reduce the likelihood of damage during removal.

5. Remove the device by simultaneously pulling the guidewire and capture sheath through the guiding catheter or interventional sheath introducer, and out of the hemostasis valve as a single unit. Care should be taken when pulling the basket through the open hemostasis valve, to avoid potential release of captured emboli.

6. Resume the interventional procedure per Preparations for Use, steps 1 through 14, and Directions for Use, steps 1 through 15.

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are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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Recommended Procedure

1. Remove guidewire from packaging, inspect and place into and/or wipe with standard heparinized saline solution.
2. Inspect and prepare selected catheter according to the manufacturer's instructions.
3. If desired to facilitate guidewire navigation through the vasculature, carefully shape the guidewire tip according to standard practices.
4. Select the catheter attachments (i.e. stopcock, hemostasis valve, etc.) as required for the indicated procedure.
5. Insert the guidewire through the catheter attachments, and/or into the guidewire port of the catheter.

Note: Free movement of the guidewire within a catheter is an essential feature of a steerable guidewire system. This movement gives the user valuable tactile information. Ensure proper guidewire size has been selected for application. Test the system for resistance prior to use.

Directions for Use

1. To facilitate vessel selection and tracking, the tip of the guidewire should be steered by carefully and slowly rotating the guidewire. For optimal torsional control of the distal guidewire tip, a torque device may be placed on the proximal portion of the guidewire. To attach the torque device, loosen the cap, slip the torque device (cap end first) over the proximal end of the guidewire, and secure it by tightening the cap firmly. Location of the torque device on the proximal guidewire section is at the discretion of the physician.
2. Confirm guidewire placement in two fluoroscopic projections to assure that the distal tip is intraluminal in the desired vessel.
3. Hold the guidewire in place while ad-

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vancing a catheter over it.

4. If a different guidewire tip configuration is desired, the distal tip of the guidewire may be reshaped according to standard practice.

Caution:

1. Attention should be focused to the guidewire movement in the vasculature. Before a guidewire is moved or torqued, the tip movement should be examined under fluoroscopy.

Do not torque a guidewire without observing corresponding movement of the tip. Always advance or withdraw a guidewire slowly. Never push, or withdraw a guidewire which meets resistance. Resistance may be detected by tactile response, and/or visual guidewire distal tip distortion observed during fluoroscopy. If resistance is determined, check the complete catheter and guidewire system.

2. Maintain continuous flush while removing and reinserting the guidewire to prevent air entry into the catheter system.

All exchanges should be done slowly to minimize air entry or trauma. When reintroducing a guidewire into a catheter within a vessel, confirm that the catheter tip is free and unobstructed within the vessel lumen.

Do not reuse. Discard after one procedure.

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Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused. Accordingly, Cordis Corporation will not be responsible for any direct, incidental or consequential damages resulting from reuse of the product.

Never push, auger, withdraw, or torque a guidewire that meets resistance. **STOP** the procedure. **DO NOT** move or torque the guidewire. First, using fluoroscopy, determine the cause of resistance and take any necessary remedial action. Torquing or pushing a guidewire against resistance may cause guidewire damage, and/or guidewire tip separation, or direct damage to the vessel. Resistance may be felt and/or observed (via fluoroscopy) by noting any buckling of the guidewire tip. If guidewire tip prolapse is observed, **DO NOT** allow the tip to remain in a prolapsed position; otherwise damage to the guidewire may occur.

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Instruction for Use Guidewires FrontRunner XP

The FRONTRUNNER XP CTO consists of a handle assembly with an integral rotator, a flush port, a flexible shaft, and a radiopaque blunt-shaped distal actuating tip. The handle lever provides manual opening and closing of the distal tip, and the handle rotator provides rotational control for the shaft and distal tip.

The FRONTRUNNER XP CTO is intended to facilitate the intraluminal placement of conventional guide wires beyond stenotic lesions (including chronic total occlusions) in the peripheral vasculature prior to further percutaneous intervention

This device is not intended for use in the cerebral vasculature.

This product is designed for use by physicians trained in and familiar with percutaneous interventional techniques. Possible complications may include, but are not limited to, the following:

- vessel dissection, perforation, or injury
- embolism
- vascular thrombosis
- pseudoaneurysm
- vessel spasm • infection

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- Do not use without completely reading and understanding this document.
- Store in a cool, dark, dry place.
- Do not expose to organic solvents (e.g., alcohol).
- Excessive bending or kinking may affect performance.
- Torquing catheter excessively may cause damage to the product. Withdraw the catheter should it become kinked.
- If strong resistance is felt during manipulation, determine the cause of the resistance before proceeding further. If the cause cannot be determined, withdraw the medical device.

1. Use sterile technique to carefully remove the catheter from the packaging. Inspect the catheter to ensure the catheter exhibits no signs of damage.

2. Visually inspect for complete opening and closing of the distal tip. Push on the distal end of the handle lever to open the distal tip. Push on the proximal end of the handle lever to close the distal tip.

Do not use if full closure cannot be achieved.

3. If necessary, the distal shaft of the FRONTRUNNER XP Catheter (FR-XP) may be manually shaped (curved) to facilitate its guidance through the vasculature. Do not use tools of any kind to shape the FR-XP as this may damage the FR-XP.

4. Flush the FR-XP with heparinized saline through the proximal flush port located on the handle assembly.

5. Wipe the FR-XP with heparinized saline or submerge in heparinized saline to hydrate the hydrophilic coating just prior to use.

6. Prepare the Micro Guide Catheter XP (MGC) according to its Instructions for Use.

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7. Advance the MGC over the guide wire to the lesion. Remove the guide wire.

8. Introduce the FR-XP through the MGC.

Alternatively, the FR-XP may be loaded into the MGC prior to introducing the FR-XP/MGC system into the vasculature.

9. With the distal tip in the closed position, advance the FR-XP to the lesion using fluoroscopic guidance and steering the FR-XP distal tip with the handle rotator.

10. Engage distal tip in lesion: Continue to advance the closed distal tip into the lesion, then open the distal tip by pushing on the distal end of the handle lever.

11. Disengage distal tip from lesion: Close the distal tip by pushing on the proximal end of the handle and retract the FR-XP approximately 1 cm. Do not hold the handle closed, or torque the catheter, while retracting it. Ensure that the tip is disengaged from tissue. If the FR-XP does not retract freely, continue to retract without torquing the FR-XP. Once disengaged, move the handle to the fully closed position. Using high-quality fluoroscopy ensure that the distal tip is fully closed.

Do not rotate or advance the FR-XP until disengagement from the lesion and full closure of the distal tip is confirmed by fluoroscopy.

12. Re-engage the distal tip in the lesion: Repeat steps 9, 10 and 11 and advance the FR-XP as needed, to create a path through the occlusion. Reorient the closed distal tip, as required, by turning the rotator on the handle assembly. Monitor all distal tip opening and closing using fluoroscopic guidance.

13. Upon crossing the occlusion, advance the MGC over the FR-XP past the occlusion. Use fluoroscopy to ensure the distal tip is fully closed. Carefully retract the FR-XP completely.


14. Introduce a guide wire and advance it through the MGC. Proceed with further percutaneous intervention.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no cir-

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cumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

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Catheter Sheath Introducer

A catheter sheath introducer (CSI) facilitates percutaneous entry of an intravascular device. Color coding below indicates the French size of the largest intravascular device that will pass through the CSI. The availability of French sizes varies per product line. Refer to product catalog for availability.

Note: Color coding not assigned to 12F, 13F, and 14F catheter sheath introducers.

French Color

- 4 (1.35 mm) Red
- 5 (1.65 mm) Gray 6 (2.0 mm) Green 7 (2.3 mm) Orange
- 8 (2.7 mm) Blue 9 (3.0 mm) Black
- 10 (3.3 mm) Magenta
- 11 (3.7 mm) Yellow
- 12 (4.0 mm) N/A
- 13 (4.3 mm) N/A
- 14 (4.7 mm) N/A

Vessel Dilator (Included or available separately)

A vessel dilator facilitates the percutaneous entry of the CSI by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

Mini-Guidewire (Included or available separately)

A mini-guidewire is provided for maintaining access to the puncture site upon removal of the needle and the insertion of the CSI.

Obturator (Included or available separately)

An obturator supports the cannula wall of an indwelling sheath. Use an obturator that is one

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French size smaller than the CSI. This allows for flushing, infusing and pressure monitoring while the obturator is in place. Use chart below as a reference.

French Size

CSI 4 5 6 7 8 9 10 11

Obturator - 4 5 6 7 8 9 10

Needle (Included or available separately)

A needle provides access to the vessel and serves as a pathway for the introduction of the guidewire.

The CSI is indicated for use in arterial and venous procedures requiring percutaneous introduction of intravascular devices

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- air embolism
- infection
- intimal tear
- hematoma at the puncture site
- perforation of the vessel wall
- thrombus formation.

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For one use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired through cleaning, resterilization or reuse and may cause adverse patient reactions. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from reuse of the CSI.

This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use.

Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

Do not use with Ethiodol* or Lipiodol (Ethiodol and Lipiodol are trademarks of Guerbet S.A) contrast media, or other such contrast media, which incorporate the components of these agents.

Do not leave a CSI in place for extended periods of time without a catheter or an obturator to support the cannula wall.

- Store in a dry, dark, cool place.
- Do not use if package is open or damaged.
- Note "Use Before" or "Use By" date prior to using product.
- Do not resterilize. Exposure to tem-

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temperatures above 54°C (130°F) may damage the catheter sheath and components.

- Do not expose to organic solvents, e.g. alcohol.

- If increased resistance is felt upon insertion of the CSI, investigate the cause before continuing. If the cause of the resistance cannot be determined and corrected, discontinue the procedure and withdraw the CSI.

1. Remove CSI from package using sterile technique.

2. Remove air from CSI by flushing with heparinized saline or suitable isotonic solution.

3. Insert the vessel dilator through the CSI's hemostasis valve, snapping it into place at the hub (Figure 1). Flush with heparinized saline or suitable isotonic solution.

4. Introduce the cannula of an angiographic needle into the vessel using aseptic technique. Holding the needle in place, insert the flexible end of the guidewire through the needle and into the vessel. If a J tip is used, slide the guidewire introducer over the J to straighten it for insertion. Gently advance the guidewire to the desired depth.

Note: Refer to product labeling for the guidewire size that is compatible with the system components.

Note: If a needle with a metal cannula is used, do not withdraw the guidewire after it has been inserted as it may damage the guidewire.

5. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the CSI is inserted

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into the vasculature.

6. Thread the CSI/dilator assembly over the guidewire, grasping the sheath close to the skin to prevent buckling. Using a

rotating motion, advance the assembly through the tissue and into the vessel (Figure 2).

7. Detach the vessel dilator from the CSI by releasing the snap-fit ring at the hub

(Figure 3). Withdraw the guidewire and dilator. To avoid damage to the CSI tip, do not withdraw the dilator until the CSI is in the vessel.

8. Aspirate from the sideport extension to remove any potential air. After aspiration

and flushing, consider establishing a heparinized solution or suitable isotonic solution via the sideport extension. The addition of a heparinized saline drip via the sideport can help in the prevention of thrombus formation.

9. Introduce the selected catheter into the CSI using one of the following methods:

а. Straighten the catheter tip by hand.

б. Insert a guidewire into the catheter until the tip is straight.

Note: Hold the sheath in place when inserting, positioning, or removing the catheter. The suture collar (Figure 1) may be used as a temporary suture site.

Note: If there is a reduction in the valve hemostasis, insert the tip of the vessel dilator into the valve and withdraw it slowly.

Note: When measuring right atrial pressure or determining cardiac output by thermodilution methods, temporarily discontinue infusion through the sideport extension to help prevent errors in hemodynamic measurement.

10. To exchange catheters, slowly withdraw the catheter from the CSI and repeat the insertion procedure.

11. Remove the sheath when clinically indicated by placing compression on the vessel above the puncture site, and slowly withdraw the CSI. Discard the

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sheath appropriately.

Important: Upon removing any CSI, aspirate via the sideport extension to collect any fibrin that may have been deposited within or at the tip of the sheath.

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
1. Flush the CSI with heparinized saline. Holding the CSI in place, slide the obturator through the hemostasis valve pushing the snap-fit ring into place. A heparinized solution or suitable isotonic solution infused through the sideport extension is recommended.

Note: Use an obturator that is one French size smaller than the CSI to allow for flushing, infusion, and pressure monitoring. Refer to table in “Obturator” section under System Components / Description.

2. To remove, hold the CSI in place, releasing the snap-fit ring, and withdraw.

3. Discard the sheath and obturator accordingly.

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Contents (Included or available separately)

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After use all used components and packing materials may be potentially biologically dangerous. Their treatment and disposal are performed according to the accepted medical practice, and also according to local, state and federal laws and provisions. The products which had a contact with blood and / or body fluids are subject to disposal and / or destruction ac-

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according to SanPiN 2.1.7.2790-10 (B class). Unused products (which had no contact with blood and / or body fluids), including those with expired shelf-life, are subject to disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

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The Cordis Long Sheath is an assembly consisting of:

- A radiopaque polyethylene sheath
- A 22-cm sideport extension with stopcock for flushing and aspirating the sheath
- A hemostasis valve* for minimizing back-bleeding.

The Cordis Long Sheath Set consists of the Long Sheath, a catheter, and a vessel dilator. All components that enter blood vessels are radiopaque.

The Long Sheath and Long Sheath Set are available in both straight and curved configurations and in two lengths. The 43-cm sheath (recommended for use with 50-cm Cordis Biopsy Forceps) may be used in the jugular vein for right ventricular biopsies, or to introduce catheters into iliac arteries. The 96-cm sheath (recommended for use with 104-cm Cordis Biopsy Forceps) may be used in the femoral artery or vein for right or left ventricular biopsies. The color coding below indicates the French size of the Cordis catheter that should be used for introduction. The catheter is inserted through the sheath, forming an assembly to be introduced as a unit using the Seldinger technique.

French Color

7	Orange
8	Blue

Catheters, guidewires, and angiographic needles are available from Cordis.

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to provide percutaneous access of intravascular devices. The Long Sheath Set are indicated for use in arterial and venous procedures requiring percutaneous introduction of Cordis Biopsy Forceps and catheters.

None known.

Possible complications include air embolism, perforation of the vessel wall, infection, and hematoma at the puncture site.

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Discard after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. All parts are extremely difficult to clean adequately after exposure to biological materials and may cause adverse patient reactions if reused. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from reuse of the Long Sheath or Long Sheath Set.

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Caution: To reduce the possibility of air aspiration through the sheath during intravenous procedures, advise the patient to maintain an appropriate position and to avoid sudden movements.

1. Introduce the cannula of an angiographic needle into the vessel (Figure 1).

Note: See the Cordis Angiographic Product Systems Pricing Catalog for the correct needle size.

2. Holding the cannula in place, insert the flexible end of a guidewire approximately 15 cm into the vessel (Figure 2).

Note: Guidewire diameter must match the lumen of the catheter. See the catheter hub or package label for the appropriate guidewire diameter.

Note: If a guidewire is used with a needle having a metal cannula, do not withdraw the guidewire after it has been inserted. Withdrawal may damage the guidewire.

3. Hold the guidewire in place and withdraw the needle.

4. Slide a vessel dilator over the guidewire, into the tissue, and through the vessel wall to create a pathway for the catheter tip (Figure 3).

Note: The vessel dilator must be one French size larger than the catheter.

5. Holding the guidewire in place, slide the dilator off the guidewire and discard.

6. Fill a syringe with heparinized saline solution. (If heparin is contraindicated, use another isotonic solution.) Attach the sy-

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ринге to the sideport extension stopcock.

7. Invert the sheath assembly and flush with solution to remove air.

8. Insert the catheter into the sheath assembly (Figure 4) and flush the sheath and catheter with the solution or with contrast medium.

9. Remove the syringe from the sideport extension stopcock and replace it with a saline solution IV drip.

10. Thread the catheter/sheath assembly over the guidewire and into the vessel, being careful that the catheter does not slide back through the sheath (Figure 5).

Important: Maintain an IV drip whenever inserting or manipulating a catheter to prevent backflow of blood into the sheath.

Note: A guidewire may be used to navigate tortuous vessels or to position the catheter tip for selective angiographic procedures. If an obstruction is encountered, contrast medium may be injected through the catheter to help visualize the area without removing the catheter or sheath assembly.

11. Withdraw the guidewire from the catheter/sheath assembly.

Note: If strong resistance is felt when removing the guidewire, remove the catheter and guidewire as a unit to prevent injury to the blood vessel and damage to the catheter and guidewire.

12. Position the catheter/sheath assembly in the desired location. At this time, contrast medium may be injected through the catheter to opacify the vessels.

Note: Hold the sheath assembly in place whenever inserting, positioning, or removing a catheter. The collar on the sideport extension may be used as a temporary suture site.

13. Remove the catheter from the assembly.

Important: When removing the catheter, aspirate the sideport extension with a syringe to collect any fibrin that may have been deposited at the tip of the sheath.

Do not aspirate the sideport extension after removing the catheter. If aspiration becomes necessary, tightly cover the hemostasis valve with the valve cap or with a

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


thumb to help prevent air intrusion through the valve.

The sheath assembly can now be used for the introduction of Cordis Biopsy Forceps or another catheter.

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Instruction for Use Biopsy forceps

The radiopaque Cordis Disposable Biopsy Forceps consist of a 3-pull ring handle, stainless steel cutting jaws, and a coiled shaft. The flexible thumb ring rotates to accommodate any position.

A spring in the handle keeps the cutting jaws closed. The jaws are opened by moving the double rings away from the thumb ring.

The jaws are closed on the specimen by moving the double rings toward the thumb ring. The forceps can remove a tissue sample of approximately 2.46 mm³ for the 5.5F (1.85 mm) and 5.20 mm³ for the 7F (2.3 mm) catalog items.

The 50 cm biopsy forceps are intended for right ventricular biopsies using the internal jugular (IJ) and subclavian (SC) approach.

The 104 cm biopsy forceps are intended for right or left ventricular biopsies using the femoral (FEM) approach.

A Cordis 6F (2.0 mm) sheath introducer is recommended for the introduction of the 5.5F biopsy forceps. A Cordis 7F sheath introducer is recommended for the introduction of the 7F biopsy forceps.

The Cordis biopsy forceps include color coding of the handle strain relief.

The color coding system is as follows:

6F (2.0 mm) Green

7F (2.3 mm) Orange

The Cordis biopsy forceps are designed for endomyocardial biopsies.

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- hematoma at the puncture site
- infection
- perforation of the vessel wall or the myocardium
- vessel trauma
- embolism
- death.

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Discard forceps after completing one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Forceps are extremely difficult to clean after exposure to biological material and may cause adverse patient reactions if reused.

- The forceps should be thoroughly rinsed with heparinized saline before and after each biopsy during the procedure.
- The heart should be routinely monitored by ECG during the procedure.
- Use prior to the "Use By" date.

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Note: These biopsy forceps are NOT formable.

1. Refer to the "Description" section for aid in selecting the size of the sheath introducer and biopsy forceps to be used.

2. Using fluoroscopy, advance the forceps through the appropriate sheath introducer and into the right or left ventricle.

Note: When using the 104 cm biopsy forceps, confirm the position of the sheath introducer in the ventricle prior to each biopsy to be performed.

3. After confirming that the tip of the forceps is in the ventricle, open the jaws.

4. Advance the open jaws to the heart wall. Close the jaws firmly to obtain a tissue specimen. Maintain sufficient pressure on the double rings to assure retention of specimen during withdrawal.

Note: Closing of the jaws and withdrawal of the bioptome should be performed in a single motion.

5. Flush the sheath introducer continuously with heparinized saline while slowly withdrawing the forceps.

6. Remove the sample from jaws.

7. Rinse bioptome before and after each biopsy.

8. Repeat steps 2 thru 7 until sufficient tissue volume is obtained.

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The radiopaque Cordis Formable/Torqueable Disposable Biopsy Forceps consist of a 3-pull ring handle, stainless steel cutting jaws, a coiled shaft, and PTFE encased formable shaft. The flexible thumb ring rotates to accommodate any position. The forceps (not including handle) have a diameter of 2.3 mm.

A spring in the handle keeps the cutting jaws closed. The jaws are opened by moving the double rings away from the thumb ring. The jaws are closed on the specimen by moving the double rings toward the thumb ring. The forceps can remove a tissue sample of approximately 5.03 mm³.

The 50 cm biopsy forceps are intended for right ventricular biopsies using the internal jugular (IJ) or subclavian (SC) approach.

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A Cordis 7F (2.3 mm) sheath introducer is recommended for the introduction of the 7F biopsy forceps.

The Cordis biopsy forceps include color coding of the handle strain relief. The color coding system is as follows:

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Procedures requiring biopsy forceps should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure.

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- hematoma at the puncture site
- infection
- perforation of the vessel wall or the myocardium
- vessel trauma
- embolism
- death

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1. Refer to the "Description" section for aid in selecting the size of the sheath introducer and biopsy forceps to be used.
2. Using fluoroscopy, advance the forceps through the appropriate sheath introducer and into the right or left ventricle.

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Note: When using the 104 cm biopsy forceps, confirm the position of the sheath introducer in the ventricle prior to each

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biopsy to be performed.

3. After confirming that the tip of the forceps is in the ventricle, open the jaws.
4. Advance the open jaws to the heart wall. Close the jaws firmly to obtain a tissue specimen. Maintain sufficient pressure on the double rings to assure retention of specimen during withdrawal.
Note: Closing of the jaws and withdrawal of the bioptome should be performed in a single motion.
5. Flush the sheath introducer continuously with heparinized saline while slowly withdrawing the forceps.
6. Remove the sample from jaws.
7. Rinse bioptome before and after each biopsy.
8. Repeat steps 2 thru 7 until sufficient tissue volume is obtained.

Note: These biopsy forceps are formable at the tip. To modify the angle of the existing curve, grasp the distal end between the thumb and index finger (Figure 1). Slide fingers along tip while bending to desired angle (Figure 2 and Figure 3).

Discard the forceps after completing one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Forceps are extremely difficult to clean after exposure to biological material and may cause adverse patient reactions if reused.

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I, **Miguel Ávila**, in my capacity as Vice President, Global Quality and Regulatory Affairs – Cordis Corporation certify that the attached **Instruction for Use for Vessel Dilators /Obturators** issued by **Cordis Corporation** is true and correct to the best of my knowledge.

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On April 20, 2017, before me, ZHENG WANG, Notary Public, personally appeared Miguel Ávila, who proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument and acknowledged to me that she executed the same in her authorized capacity, and that by her signature on the instrument the person, or the entity upon behalf of which the person acted, executed the instrument.

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Instruction for Use Dilators/Obturator

Vessel Dilator

A vessel dilator facilitates the percutaneous entry of the catheter sheath introducer (CSI) by forming an atraumatic transition from the skin through the subcutaneous tissue to the vessel.

Note: The vessel dilator size must match the guidewire, sheath, and catheter. Refer to the product labeling or hub for size information.

Obturator

An obturator supports the cannula wall of an indwelling sheath. Use an obturator that is one French size smaller than the CSI. This allows for flushing, infusing and pressure monitoring while the obturator is in place. Use chart below as a reference.

French Size

CSI 4 5 6 7 8 9 10 11

Obturator - 4 5 6 7 8 9 10

Color coding below refers to the French size of the largest intravascular device that will pass through the CSI. Color-coding of the vessel dilator matches that of the corresponding CSI. Color-coding will indicate the actual French size of the obturator which is one French size smaller than the CSI.

The availability of French sizes varies per product line. Refer to product catalog for availability.

French	Color	French	Color
4 (1.35 mm)	Red	8 (2.7 mm)	Blue
5 (1.65 mm)	Gray	9 (3.0 mm)	Black
6 (2.0 mm)	Green	10 (3.3 mm)	Magenta
7 (2.3 mm)	Orange	11 (3.7 mm)	Yellow

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The obturator and vessel dilator are used in conjunction with the Cordis catheter sheath introducer system, indicated for use in arterial and venous procedures requiring percutaneous introduction of intravascular devices.

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Possible complications include, but are not limited to:

- air embolism
- infection
- intimal tear
- hematoma
- perforation of the vessel wall
- thrombus formation

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- For one use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired through cleaning, resterilization or reuse and may cause adverse patient reactions.
- Do not use with Ethiodol or Lipiodol* contrast media, or other such contrast media, which incorporate the components of these agents.
- Do not leave a CSI in place for extended periods of time without a catheter or obturator to support the cannula wall.

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Vessel Dilator Insertion and Withdrawal

1. To dilate, slide the vessel dilator over the guidewire and into the tissue and vessel wall to create a pathway for the catheter tip.
2. While holding the guidewire in place, gently withdraw the dilator from the entry site and off the guidewire. Discard appropriately.

Obturator Insertion and Withdrawal

1. Flush the CSI with heparinized saline. Holding the CSI in place, slide the obturator through the hemostasis valve pushing the snap-fit ring into place. A heparin drip through the sideport extension is recommended. If heparin is contraindicated, use another suitable isotonic solution.
2. To remove, hold the CSI in place, releasing the snap fit ring, and withdraw. Discard appropriately.

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Instruction for Use Stents PRECISE PRO RX

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
The Cordis PRECISE PRO Rx Nitinol Stent System consists of a nitinol self-expanding stent preloaded on a .065" (1.65 mm) or .078" (1.98 mm) sheathed delivery system. The delivery system consists mainly of an inner shaft and an outer sheath with radio-opaque markers, and a Tuohy Borst valve. The inner shaft consists of a support member and wire lumen. The proximal portion of the support member is comprised of a hub connected to a stainless steel wire and hypotube and distally of a stainless steel coil. The wire lumen originates distally in a catheter tip and terminates proximally at a guidewire exit port designed to accept a .014" (0.36 mm) guidewire. The outer sheath has a proximal shaft and distal outer sheath with a nominal working length of 135 cm. The self-expanding PRECISE stent is constrained within the space between the inner shaft and the distal outer sheath, located between distal and proximal stent markers on the inner shaft. The stent expands to its unconstrained diameter when released from the deployment catheter into the vessel. Upon deployment, the stent forms an open lattice and pushes outward on the luminal surface, helping to maintain the patency of the vessel. Due to the self-expanding behavior of nitinol, the stents are indicated for placement into vessels that are 1–2 mm smaller in diameter than the unconstrained diameter of the stent. Device depictions and components are provided in Figure 1.

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The Cordis **PRECISE PRO Rx** Nitinol Stent System is indicated for use in patients with stenotic lesions of the carotid artery(ies).

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Generally, contraindications to PTA are also contraindications for stent placement. Contraindications include, but may not be limited to:

- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders.
- Stenting of intra-cranial arteries.

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Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications.

Complications may occur at any time during or after the procedure.

Potential complications may include, but are not limited to:

- Death
- Respiratory arrest
- Emergency artery bypass graft surgery
- Hemorrhagic or embolic stroke / TIA
- Renal failure
- Sepsis / infection
- Embolism
- Coronary ischemia
- Arrhythmia

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- Drug reactions, allergic reaction to contrast medium or to the implanted device
- Vascular injury, including perforation, rupture, and dissection
- Disseminated intravascular coagulation
- New or worse encephalopathy
- G.I. bleeding from anticoagulation / antiplatelet medication
- Hemorrhage
- Parenchymal hemorrhage
- Aneurysm and pseudoaneurysm formation
- Intimal tear / dissection
- Stent migration / embolization
- Thrombosis
- Bradycardia and hypotension
- Arteriovenous fistula
- Tissue necrosis
- Stent misplacement
- Vessel occlusion, restenosis or recurrent stricture
- Hematoma
- Carotid artery spasm

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Pre- Procedure

The patient should be started on nonbuffered, nonenteric-coated aspirin 72 hours prior to the procedure per standard hospital dosing guidelines or as prescribed by a physician. Antiplatelet therapy should be administered 24-48 hours

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prior to the procedure, according to hospital protocol. Antiplatelet therapy following a carotid stenting procedure should be administered per physician instructions.

The percutaneous placement of the stent in a stenotic carotid artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. Access vessels must be sufficiently patent or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

1. Inject Contrast Media
Perform a percutaneous angiogram using standard technique.
2. Identify and mark the lesion
Fluoroscopically identify and mark the lesion, observing the most distal level of the stenosis.

Device Selection and Preparation

1. Select Stent Size Measure the length of the target lesion to determine the length of stent(s) required. Measure the diameter of the reference vessel (proximal and distal to the lesion). It is necessary to select a stent which has an unconstrained diameter that is at least 1 mm larger than the largest reference vessel diameter to achieve secure placement according to the following Stent Size Selection Table.

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Stent Size Selection Table

Vessel Lumen Diameter	Unconstrained Stent Diameter	% Length Foreshortening
3.0–4.0 mm	5.0 mm	1.2%
4.0–5.0 mm	6.0 mm	2.4%
5.0–6.0 mm	7.0 mm	4.1%
6.0–7.0 mm	8.0 mm	6.2%
7.0–8.0 mm	9.0 mm	5.8%
8.0–9.0 mm	10.0 mm	8.0%

Refer to product labeling for stent length.

Note: The percent foreshortening of stent length is based upon a mathematical calculation.

2. Preparation of Stent Delivery System

CAUTION: The stent delivery system is shipped with the Tuohy Borst valve OPEN. Be careful not to prematurely deploy the stent during preparation. Prep the device in the tray per instructions below. Close the Tuohy Borst valve prior to removing the device from the tray.

a. Open the outer box to reveal the pouch containing the stent and delivery system.

b. Check the temperature exposure indicator on the pouch to confirm that the black dotted pattern with a grey background is clearly visible. See Warnings section.

c. After careful inspection of the pouch looking for damage to the sterile barrier, carefully peel open the pouch and remove the tray. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

d. While in the tray, attach a stopcock to the Y connection on the Tuohy Borst valve.

e. Attach a 5-cc syringe filled with heparinized saline to the open stopcock and apply positive pressure until saline weeps from the proximal end of the Tuohy Borst valve. Lock the Tuohy Borst valve.

f. Close the stopcock attached to the Tuohy Borst Y connection.

g. Extract the stent delivery system from the tray. Examine the device for any damage. Evaluate the distal end of the catheter to ensure that the stent is contained within the outer sheath. Do not use if the stent is partially deployed. If a gap between the catheter tip and outer sheath tip exists,

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open the Tuohy Borst valve and gently pull the inner shaft in a proximal direction until the gap is closed. Lock the Tuohy Borst valve after the adjustment by rotating the proximal valve end in a clockwise direction.

Stent Deployment Procedure

1. Insertion of Introducer Sheath or Guiding Catheter and Cordis ANGIOGUARD RX Emboli Capture Guidewire System.

a. Access the treatment site utilizing the appropriately sized accessory equipment. See Table "Recommended Accessory Catheter Sizing."

Table 1: Recommended Accessory Catheter Sizing

PRECISE PRO Rx Stent Diameter	Minimum Sheath	Minimum Guiding
	Introducer	Catheter I.D.
5, 6, 7, 8 mm	5F (1.98 mm)	.078" (1.98 mm)
9 & 10 mm	6F (2.21 mm)	.087" (2.21 mm)

b. Insert an appropriately sized Cordis ANGIOGUARD RX Emboli Capture Guidewire System via the introducer sheath or guiding catheter.

NOTE: PLEASE REFER TO THE CORDIS ANGIOGUARD RX INSTRUCTIONS FOR USE FOR PLACEMENT PROCEDURE AND USE OF THE DEVICE.

c. The Cordis PRECISE PRO Rx Nitinol Stent System is compatible with a .014" (0.36 mm) or smaller guidewire.

2. Dilation of Lesion

a. If appropriate, pre-dilate the lesion using standard PTA techniques.

b. Remove the PTA balloon catheter from the patient maintaining lesion access with the guidewire.

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3. Introduction of Stent Delivery System

а. Flush the guidewire lumen of the stent delivery system with heparinized saline by connecting a 5-cc syringe filled with heparinized saline solution to the stopcock attached to the Y connection (9) on the Tuohy Borst valve (1) to expel air. Ensure that the Tuohy Borst valve (1) is in the locked position to prevent premature stent deployment. Apply positive pressure to the syringe until saline weeps from the guidewire exit port (16). While covering the guidewire exit port (16) with thumb and forefinger, apply positive pressure to the syringe until saline weeps from the catheter tip (4) and the space between the outer sheath radiopaque marker (11) and the catheter tip (4). Continue to flush to ensure all air is removed from the system, then close the stopcock attached to the Y connection (9) on the Tuohy Borst valve (1).

б. Ensure that the Tuohy Borst valve connecting the inner shaft and outer sheath is locked by rotating the proximal valve end in a clockwise direction to prevent premature stent deployment.

с. Advance the PRECISE PRO Rx System over the .014" (0.36 mm) ANGIOGUARD RX Emboli Capture Guidewire System until the guidewire exit port (16) is just outside the accessory device's Tuohy Borst valve. Adjust the accessory device's Tuohy Borst valve to maintain a snug seal about the PRECISE PRO Rx System. Look for and confirm back flow through the guidewire exit port (16) opening.

д. After confirming back flow, advance the PRECISE PRO Rx System. Again adjust the accessory device's Tuohy Borst valve to maintain a snug seal, now over the .038" (0.97 mm) proximal shaft (6A) and continue to advance the PRECISE PRO Rx System to the lesion site.

е. Prior to contrast injection, confirm again a snug seal between the .038" (0.97

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но proximal shaft (6A) of the PRECISE PRO Rx System and the accessory device's Tuohy Borst valve. FAILURE TO DO SO COULD RESULT IN AIR INTRODUCTION DURING ASPIRATION, RESULTING FROM A POOR SEAL.

NOTE: If resistance is met during delivery system introduction, the system should be withdrawn and another system should be used.

4. Slack Removal

а. Advance the stent delivery system past the lesion site.

б. Pull back the stent delivery system until the radiopaque inner shaft markers (leading and trailing ends) move in position so that they are proximal and distal to the target lesion.

с. Ensure that the stent delivery system outside the patient remains flat and straight.

CAUTION: Slack in the catheter shaft either outside or inside the patient may result in deploying the stent beyond the lesion site.

5. Stent deployment

NOTE: When ready to proceed with stent deployment, heparin may be administered per standard hospital practice or as prescribed by a physician. Heparin may be continued following stent deployment if so indicated by a physician or hospital protocol.

а. Verify that the delivery system's radiopaque inner shaft markers (leading and trailing ends) are proximal and distal to the target lesion.

б. Unlock the Tuohy Borst valve connecting the inner shaft and outer sheath of the delivery system.

с. Ensure that the access sheath or guiding catheter does not move during deployment.

д. Initiate stent deployment by retracting the outer sheath while holding the inner shaft in a fixed position. Deployment is complete when the outer sheath marker passes the proximal inner shaft stent marker.

NOTE: The mechanism for stent deployment is outer sheath

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retraction. Deployment is completed by maintaining inner shaft position while retracting the outer sheath and allowing the stent to expand (often referred to as the “pin-and-pull” method).

NOTE: When more than one stent is required to cover the lesion, the more distal stent should be placed first. Efforts should be taken to minimize the stent overlap. In no instance, should more than two (2) stents ever overlap.

6. Post-deployment Stent Dilatation

a. While using fluoroscopy, withdraw the entire delivery system as one unit, over the guidewire and out of the body. Remove the delivery device from the guidewire.

NOTE: If any resistance is met during delivery system withdrawal, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the system as one unit. (Do not remove guidewire.)

b. Using fluoroscopy, visualize the stent to verify full deployment.

c. If incomplete expansion exists within the stent at any point along the lesion, post deployment balloon dilatation (standard PTA technique) can be performed.

d. Select an appropriately sized PTA balloon catheter and dilate the lesion with conventional technique.

The inflation diameter of the PTA balloon catheter used for post dilatation should approximate the diameter of the reference vessel. Remove the PTA balloon catheter from the patient.

7. Post Stent Placement

a. A post stent angiogram should be obtained.

b. Remove the ANGIOGUARD RX Emboli Capture Guidewire System in accordance with ANGIOGUARD RX Instructions for Use. Remove the sheath and establish hemostasis.

c. Discard the delivery system, guidewire, and sheath.

NOTE: Physician experience and discretion will determine the appropriate

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post procedure drug regimen for each patient.

Through non-clinical testing, the Cordis **PRECISE** Nitinol Stent has been shown to be MRI conditional at field strengths of 3 Tesla or less and a maximum whole body averaged specific absorption rate (SAR) of 0.8W/kg for 20:00 minutes of MRI. The **PRECISE** stent should not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of stent migration at field strengths higher than 3 Tesla. In this testing, the stent produced a temperature rise of 0.5 degrees C at a maximum whole body averaged SAR of 0.8W/kg for 20:00 minutes of MRI. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent. The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.

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signed for the use of lower injection systems.

- The delivery system is not recommended for use with a leaflet type valve.
- When catheters are in the body, they should be manipulated only under fluoroscopy.
- Radiographic equipment that provides high quality images is needed.
- When treating multiple lesions, the distal lesion should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chance for dislodging stents which may have already been placed.
- Recrossing a deployed stent with adjunct devices must be performed with caution.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infection, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.
- In patients requiring the use of antacids and/or H₂-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g. aspirin) may be adversely affected.
- Venous access should be available during carotid stenting in order to manage bradycardia and/or hypotension either by pharmaceutical intervention or placing of a temporary pacemaker, if needed.
- Both PRECISE PRO Rx Nitinol Stent Systems are shipped with the hemovalve in the OPEN position (see "Preparation of the Stent Delivery System").
- Prior to stent deployment remove all slack from the catheter delivery system (see "Stent Deploy-

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ment”).

- Ensure that the catheter system is flushed according to the steps outlined in “Introduction of Stent Delivery System” and “Preparation of Stent Delivery System.” Failure to do so could result in air entering the access catheter.

- Ensure that there is a tight seal between the PRECISE PRO Rx catheter and the valve for the access catheter during aspiration. Failure to do so could result in air entering the access catheter.

- Store in a cool, dark, dry place.

The manufacturer warrants conformance of this medical device with the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer’s use instructions are met.

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Descriptions or specifications in Cordis printed matter, including this publication,

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The Cordis S.M.A.R.T. Nitinol Stent System is designed to deliver a self-expanding stent to the peripheral vasculature via a 6F (2.0 mm) sheathed delivery system. The self-expanding stent is composed of a nickel titanium alloy (nitinol). A total of 12 (6 at each end) tantalum radiopaque markers are located on the ends of the stent. The stent is a flexible, fine mesh tubular prosthesis, which expands upon deployment to appose the vessel wall. Upon deployment, the stent imparts an outward radial force on the luminal surface of the vessel to establish patency.

The delivery system, as shown in Figure 1, is comprised of an inner shaft and an outer sheath that are locked together with a Tuohy Borst valve (1). The inner shaft is comprised proximally of a stainless steel tube (2), and distally of a stainless steel coil covered with a polymeric tube (3). The inner shaft terminates distally in a catheter tip (4) and originates proximally in a Luer hub (7) designed to accept a .035" (0.89 mm) guidewire.

The 6F outer sheath (6) connects proximally to

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
the Tuohy Borst valve (1) via a Luer hub (5). The self-expanding stent (8) is constrained within the space between the inner shaft (3) and the outer sheath (6). This space is flushed prior to the interventional procedure by injecting fluid via the Y connection (9) on the Tuohy Borst valve. Stent movement during sheath retraction is restricted by a radiopaque marker (10) connected to the inner shaft. The outer sheath has a radiopaque marker (11) at its distal end.

Stent positioning about the target stricture is achieved prior to deployment utilizing the distal stent markers (12) and the proximal stent markers (13). For stent deployment, the Tuohy Borst valve is unlocked on the inner shaft by a counter clockwise rotation of the proximal valve end (14). Sheath retraction is achieved by grasping the inner shaft (7) in a fixed position and moving the outer sheath proximally relative to the inner shaft. During sheath retraction, prior to stent apposition to the wall, it may be necessary to slightly advance the entire delivery system to maintain stent positioning. Complete deployment of the stent is achieved when the proximal end of the stent and the proximal stent markers (13) visibly appose the vessel wall, and the outer sheath radiopaque marker (11) is proximal to the inner shaft stent stop (10).

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The Cordis S.M.A.R.T. Nitinol Stent System is indicated for use in patients with atherosclerotic disease of peripheral arteries, including iliac and superficial femoral, for TIPSS (Transjugular Intrahepatic Portosystemic Shunt) procedures and for palliation of malignant neoplasms in the biliary tree.

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Contraindications associated with the use of this device include:

- Bleeding disorders
- Inability to tolerate antithrombotic or antiplatelet therapies
- Intra-cranial arteries
- Known hypersensitivity to nickel titanium
- Lesions that prevents complete inflation of an angioplasty balloon or proper placement of the stent or stent delivery system
- Perforated biliary duct where leakage from the duct could be exacerbated by the stent
- Perforated vessel
- Pregnancy
- Presence of severe ascites
- Stent of a duct with total biliary occlusion which cannot be crossed by the delivery catheter
- Thrombus in or near the target lesion

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Procedures requiring percutaneous catheter introduction should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during or after the procedure. Potential complications may include, but are not limited to:

- Abrupt stent closure
- Access failure
- Allergic / anaphylactoid reaction
- Anemia
- Angina / coronary ischemia / myocardial infarction
- Aneurysm or Pseudoaneurysm in vessel or vessel access site
- Arrhythmia
- Arteriosclerosis
- Atheroembolization (Blue Toe Syndrome)
- Bile duct perforation
- Bradycardia
- Death
- Disseminated intravascular coagulation
- Edema, peripheral
- Embolization (air, plaque, thrombus, device or other)
- Fever
- Fistulazation
- Gangrene
- Gastrointestinal bleed from anticoagulation/antiplatelet medication
- Hematoma/Hemorrhage
- Hypotension/hypertension
- Infection and/or sepsis
- Intimal injury/dissection/rupture/perforation
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Limb loss

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- Liver abscess
- Pain
- Pancreatitis
- Parenchymal hemorrhage
- Renal insufficiency or failure
- Restenosis of the stented segment
- Sludge occlusion
- Stent malapposition
- Stent migration
- Stent obstruction secondary to tumor ingrowth through the stent
- Stent structure fracture
- Stroke
- Thrombosis/thrombus
- Tissue necrosis
- Tumor overgrowth at the stent ends
- Vasospasm
- Worsened claudication/rest pain

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Pre-procedure

The patient may be started on nonbuffered, nonenteric-coated aspirin one or two days prior to the procedure per standard hospital dosing guidelines or as prescribed by a physician.

The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard accepted practice. (Access vessels must be sufficiently patent, or sufficiently

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реканализ, to proceed with further intervention.) Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

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2. Identify and mark the lesion or stricture

Fluoroscopically identify and mark the lesion, observing the most distal level of the stenosis or stricture.

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Stent Size Selection Table

Vessel Lumen Diameter	Unconstrained Stent Diameter	% Length Foreshortening	% Free Area
4.0 – 5.0 mm	6.0 mm	1.2%	79%
5.0 – 6.0 mm	7.0 mm	2.0%	82%
6.0 – 7.0 mm	8.0 mm	3.1%	84%

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c. After careful inspection of the pouch looking for damage to the sterile barrier, carefully peel open the pouch and extract the stent/delivery system from the tray. Examine the device for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

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д. Flush the Tuohy Borst Y valve with heparinized saline using a 3-cc syringe to expel air. Lock the Tuohy Borst valve and continue to flush until heparinized saline weeps from the distal catheter end.

е. Flush the guidewire lumen of the stent delivery system with heparinized saline using a 10-cc syringe to expel air. Continue to flush until heparinized saline flows out of the wire lumen at the distal catheter tip.

ф. Evaluate the distal end of the catheter to ensure that the stent is contained within the outer sheath. Do not use if the stent is partially deployed. If a gap between the catheter tip and outer sheath tip exists, open the Tuohy Borst valve and gently pull the inner shaft in a proximal direction until the gap is closed. Lock the Tuohy Borst valve after the adjustment by rotating the proximal valve end in a clockwise direction.

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а. Gain access at the appropriate site utilizing the appropriate accessory equipment compatible with the 6F delivery system.

б. Insert a .035" (.89 mm) guidewire of an appropriate length through the introducer sheath or guide catheter.

2. Dilation of Lesion

а. If appropriate, pre-dilate the lesion using standard PTA techniques.

б. Remove the PTA balloon catheter from the patient maintaining lesion access with the guidewire.

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Note: No predilation is generally done with malignant biliary strictures. However, if the physician determines that predilation is necessary, standard PTA techniques may be used.

Caution: During dilation, never expand the balloon such that dissection complication could occur.

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а. Ensure that the Tuohy Borst valve, connecting the inner shaft and outer sheath, is locked by rotating the proximal valve end in a clockwise direction to prevent premature stent deployment.

б. Advance the device over the guidewire to the lesion site.

Note: If resistance is met during delivery system introduction, the system should be withdrawn and another system should be used.

Caution: Always use an appropriate size introducer sheath for the implant procedure to protect vessel and access site.

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б. Pull back the stent delivery system until the radiopaque stent markers (leading and trailing ends) move in position so that they are proximal and distal to the target lesion site.

с. Ensure the device outside the patient remains flat and straight.

Caution: Slack in the catheter shaft either outside or inside the patient may result in deploying the stent beyond the target lesion site.

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б. Ensure that the access sheath or guiding catheter does not move during deployment.

с. Unlock the Tuohy Borst valve connecting the inner shaft and outer sheath of the delivery system.

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д. Initiate stent deployment by retracting the outer sheath while holding the inner shaft in a fixed position. While using fluoroscopy, maintain position of the radiopaque stent markers relative to the targeted lesion site. Watch for the distal radiopaque markers to begin separating. Separation of the distal stent markers signals that the stent is unsheathed. Continue deploying the stent until the distal end of the stent obtains full apposition with the vessel wall. Continue deploying the stent until the proximal end of the stent obtains full apposition with the vessel wall.

Note: Failure to maintain a fixed inner shaft position or constraining the catheter shaft during deployment may result in stent compression (shortening) or elongation.

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6. Post-deployment Stent Dilatation

а. While using fluoroscopy, advance the outer sheath until the outer sheath marker contacts the catheter tip and withdraw the entire delivery system as one unit, over the guide-wire and out of the sheath introducer. Remove the delivery device from the guide-wire.

Caution: The delivery system is not designed for the use of power injection.

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с. If incomplete expansion exists within the stent at any point along the lesion, post deployment balloon dilatation (standard PTA technique) can be performed.

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catheter and dilate the lesion with conventional technique. The inflation diameter of the PTA balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the PTA balloon from the patient.

7. Post Stent Placement

а. A post stent angiogram should be obtained.

б. Remove the guidewire and sheath and establish hemostasis.

Note: Physician experience and discretion will determine the appropriate post-procedure drug regimen for each patient.

Non-clinical testing has demonstrated that the S.M.A.R.T. Stent is MR Conditional for a single implanted stent.

In non-clinical testing, the stent produced a temperate rise dependent upon total stent length and magnetic field strength in a Philips Medical Systems Intera MR scanner for maximum whole-body-average specific absorption rate (SAR) of 2.0 W/kg for 15–20 minutes of scanning. The effects of localized MRI-induced heating in the clinical environment are not known.

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the S.M.A.R.T. Stent. It may be necessary to optimize MR imaging parameters for the presence of this metallic implant.

The S.M.A.R.T. Stent should not migrate in a 3 Tesla, or less, MRI environment. Non-clinical testing of two overlapping 8 x 150 mm S.M.A.R.T. Stents in a 3 Tesla Philips Medical Systems Intera MRI Scanner demonstrated magnetically induced torque and deflection force on the stents at levels less than those im-

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posed by a patient's daily activity. With regard to stent migration and torque, MRI at 3 Tesla or less may be performed immediately following the implantation of the stent. Non-clinical testing has not been performed to rule out the possibility of stent migration at field strengths above 3 Tesla.

- Appropriate antiplatelet/anticoagulant therapy should be administered pre and post procedure.
- The long term outcome following repeat dilatation of endothelialized stents is unknown at present.
- Persons allergic to nitinol (nickel titanium) may suffer an allergic reaction to this implant.
- Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- The black dotted pattern on the grey temperature exposure indicator, found on the inner pouch, must be clearly visible. Do not use if entire temperature exposure indicator is completely black as the unconstrained stent diameter may have been compromised.
- Stenting across a major bile duct branch could lead to compromised future diagnostic or therapeutic procedures.
- Any secondary expansion of the biliary stent should be performed with a device indicated for biliary stent placement.
- Do not use with Ethiodol or Lipiodol (Ethiodol and Lipiodol are trademarks of Guerbet S.A.) contrast media.

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This product should only be used by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for interventional procedures should be employed.

- Do not use power injection systems with the delivery system.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately pre-medicated.
- Use of this device requires fluoroscopic control.
- The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.
- Safety and effectiveness has not been demonstrated in:
 - o Lesions that are either totally or densely calcified
 - o Pediatric patients

System Handling - Precautions

- Store at ambient room conditions out of direct sunlight.
- Do not use beyond the "Use By" date.
- Carefully inspect the sterile packaging and device prior to use. Do not use if it appears damaged.
- Do not expose the delivery system to organic solvents.
- Do not attempt to remove the stent from the delivery system before use.
- Do not use if the outer sheath of the device cannot be flushed.
- Do not use if the inner lumen of the guide-wire tube cannot be flushed.

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Stent Placement - Precautions

• If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to the stent or lumen.

• If resistance occurs during movement through the sheath, carefully withdraw the stent system.

• If resistance is felt when initially retracting the outer deployment sheath, do not force deployment. Carefully withdraw the stent system without deploying the stent.

• Once stent deployment has begun, the stent must be fully deployed. The system is not designed for stent repositioning or recapturing.

• Use caution when crossing a deployed stent with any adjunct device.

• Prior to stent deployment, remove all slack from the catheter delivery system.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind

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Cordis to any representation or warranty except as specifically set forth herein.

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The Cordis S.M.A.R.T. CONTROL Nitinol Stent System is designed to deliver a self-expanding stent to the peripheral vasculature via a sheathed delivery system. The outer diameter of the delivery system is specified in

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the following Stent Delivery System (SDS) Table.

Stent Delivery System (SDS)	
Unconstrained Stent Diameter	SDS Outer Diameter
6–10 mm	6F (2.0 mm)
12–14 mm	7F (2.3 mm)

Refer to product labeling for stent diameter and SDS outer diameter.

The self-expanding stent is composed of a nickel titanium alloy (nitinol). A total of 12 (6 at each end) tantalum radiopaque markers are located on the ends of the stent. The stent is a flexible, fine mesh tubular prosthesis, that expands upon deployment to appose the vessel wall. Upon deployment, the stent imparts an outward radial force on the luminal surface of the vessel to establish patency.

The delivery system, as shown in Figure 1, is comprised of an inner shaft and an outer sheath that are locked together with a flushing valve (1). The inner shaft is comprised of a polymeric tube (2) covered proximally by a stainless steel hypotube (3) and, distally by a stainless steel coil (4). The inner shaft terminates distally in a catheter tip (5) and originates proximally in a Luer hub (6) designed to accept a .035" (0.89 mm) guidewire.

The outer sheath (7) connects proximally to the flushing valve (1) via a second Luer hub (8). The self-expanding stent (9) is constrained within the space between the polymeric tube (2) and the outer sheath (7). This space is flushed prior to the interventional procedure by injecting fluid via the flushing valve (1). Stent movement during sheath retraction is restricted by an inner shaft stent stop (10) connected to the inner shaft.

The outer sheath has a radiopaque marker band (11) at its distal end. Stent positioning about the target lesion is achieved prior to deployment utilizing the distal stent markers (12) and the proximal stent markers (13). For stent deployment, the locking pin (14) must be removed. Sheath retraction is achieved by grasping the handle (15) in a fixed position with the tuning dial (16) held between the thumb and index finger. Deployment is initiated by rotating the tuning dial (16) with the

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thumb and index finger [Figure 2] in a clockwise direction until the distal stent markers (12) and the distal end of the stent, visibly appose the vessel wall. With the distal stent markers (12) and the distal end of the stent apposing the vessel wall, stent deployment continues by pulling back on the deployment lever (17) [Figure 3]. Complete deployment of the stent is achieved when the proximal end of the stent and the proximal stent markers (13), visibly appose the vessel wall and the outer sheath radiopaque marker (11) is proximal to the inner shaft stent stop (10).

The Cordis S.M.A.R.T. CONTROL Nitinol Stent System is indicated for use in patients with atherosclerotic disease of peripheral arteries, including iliac and superficial femoral, for TIPSS (Transjugular Intrahepatic Portosystemic Shunt) procedures and for palliation of malignant neoplasms in the biliary tree.

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- Abrupt stent closure
- Access failure
- Allergic / anaphylactoid reaction
- Anemia
- Angina / coronary ischemia / myocardial infarction
- Aneurysm or Pseudoaneurysm in vessel or vessel access site
- Arrhythmia
- Arteriosclerosis
- Atheroembolization (Blue Toe Syndrome)
- Bile duct perforation
- Bradycardia
- Death
- Disseminated intravascular coagulation
- Edema, peripheral
- Embolization (air, plaque, thrombus, device or other)

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- Fever
- Fistulazation
- Gangrene
- Gastrointestinal bleed from anticoagulation/antiplatelet medication
- Hematoma/Hemorrhage
- Hypotension/hypertension
- Infection and/or sepsis
- Intimal injury/dissection/rupture/perforation
- Ischemia requiring intervention (bypass or amputation of toe, foot or leg)
- Limb loss
- Liver abscess
- Pain
- Pancreatitis
- Parenchymal hemorrhage
- Renal insufficiency or failure
- Restenosis of the stented segment
- Sludge occlusion
- Stent malapposition
- Stent migration
- Stent obstruction secondary to tumor ingrowth through the stent
- Stent structure fracture
- Stroke
- Thrombosis/thrombus
- Tissue necrosis
- Tumor overgrowth at the stent ends
- Vasospasm
- Worsened claudication/rest pain

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Pre-Procedure

The patient may be started on nonbuffered,

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nonenteric-coated aspirin one or two days prior to the procedure per standard hospital dosing guidelines or as prescribed by a physician.

The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard accepted practice. (Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention.) Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

Procedure

1. Inject Contrast Media

Perform a percutaneous angiogram using standard technique. For biliary procedures the biliary tree may be injected.

2. Evaluate and Mark Lesion or Stricture

Fluoroscopically evaluate and mark the lesion or stricture, observing the most distal level of the stenosis or stricture.

3. Select Stent Size

Measure the length of the target lesion to determine the length of stent(s) required.

Measure the diameter of the reference vessel (proximal and distal to the lesion). It is necessary to select a stent that has an unconstrained diameter at least 1 mm larger than the largest reference vessel diameter to achieve secure placement according to the following Stent Size Selection Table.

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Refer to product labeling for stent diameter and length. Note: The percent foreshortening of stent length is based upon mathematical calculation.

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4. Preparation of Stent Delivery System

a. Open the box to reveal the pouch containing the stent and delivery system.

b. Check the temperature exposure indicator on the pouch to confirm that the black dotted pattern with a grey background is clearly visible. See “Warnings” section.

c. After careful inspection of the pouch looking for damage to the sterile barrier, carefully peel open the pouch and extract the stent delivery system from the tray. Examine the device for any damage. If it is suspected that the sterility or performance of the device has been compromised, the device should not be used.

d. Flush the flushing valve of the stent delivery system with heparinized saline using a 3-cc syringe to expel air. Continue to flush until heparinized saline weeps from the distal catheter end.

e. Flush the guidewire lumen of the stent delivery system with heparinized saline using a 20-cc syringe to expel air. Continue to flush until heparinized saline flows out of the wire lumen at the distal catheter tip.

f. Evaluate the distal end of the catheter to ensure that the stent is contained within the outer sheath. Do not use if the stent is partially deployed.

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Stent Deployment Procedure

1. Insertion of Introducer Sheath or Guiding Catheter and Guidewire

a. Gain access at the appropriate site utilizing the appropriate accessory equipment compatible with the stent delivery system.

b. Insert an .035" (0.89 mm) guidewire of an appropriate length through the introducer sheath or guide catheter.

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Note: No predilation is generally done with malignant biliary strictures. However, if the physician determines that predilation is necessary, standard PTA techniques may be used. Caution: During dilation, never expand the balloon such that dissection complication could occur.

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a. Ensure locking pin is still in place.

Note: If the locking pin is not in place, system performance may be compromised and another system should be used.

b. Advance the stent delivery system over the guidewire through the sheath introducer to the lesion site.

Note: If resistance is met during delivery system introduction, the system should be withdrawn and another system should be used. Caution: Always use an appropriate size introducer sheath for the implant procedure to protect vessel and access site.

4. Slack Removal

a. Advance the stent delivery system past the lesion site.

b. Pull back the stent delivery system until the radiopaque stent markers (leading and trailing ends) move in position so that

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they are proximal and distal to the lesion site.

с. Ensure that the stent delivery system outside the patient remains flat and straight. Caution: Slack in the catheter shaft either outside or inside the patient may result in deploying the stent beyond the lesion site.

5. Stent Deployment

а. Verify that the delivery system's radiopaque stent markers (leading and trailing ends) are proximal and distal to the target lesion.

б. Ensure that the introducer sheath does not move during deployment.

с. Remove locking pin from handle.

д. Initiate one-handed stent deployment by rotating the tuning dial with thumb and index fingers in a clockwise direction (direction of arrow) while holding the handle in a fixed position [Figure 2].

е. While using fluoroscopy, maintain position of the radiopaque stent markers relative to the targeted lesion site. Watch for the distal radiopaque markers to begin separating. Separation of the distal stent markers signals that the stent is unsheathed. Continue turning the tuning dial until the distal end of the stent obtains full apposition with the vessel wall.

Note: Only the initial 40 mm of the stent may be unsheathed using the tuning dial.

ф. With maintaining a fixed handle position, pull back the deployment lever to unsheath the remainder of the stent [Figure 3]. Note: Failure to maintain a fixed handle position or constraining the catheter shaft during deployment may result in stent compression (shortening) or elongation.

Note: The stent may be deployed using two hands ("Pin and pull" method) by holding the proximal end of the handle stationary with one hand and sliding the deployment lever back towards the stationary hand [Figure 4].

г. Stent deployment is complete when the proximal markers appose to the vessel wall and the outer sheath radiopaque marker is proximal to the inner shaft stent stop.

Note: When more than one stent is required to cover the lesion, the more distal stent should be placed first. Efforts should be taken to minimize the stent overlap.

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6. Post-deployment Stent Dilatation

а. Advance the deployment lever to its pre-deployment position [Figure 1] while maintaining the handle in a fixed position. Recover the delivery system by pushing the lever as far forward as possible and then by turning the dial counter-clockwise, while keeping pressure on the lever, until the lever reaches the end of the slot and the tip is re-sheathed. While using fluoroscopy, withdraw the entire delivery system as one unit, over the guidewire and out of the sheath introducer. Remove the delivery device from the guidewire. Caution: The delivery system is not designed for the use of power injection.

б. Using fluoroscopy, visualize the stent to verify full deployment.

в. If incomplete expansion exists within the stent at any point along the lesion, post-deployment balloon dilatation (standard PTA technique) can be performed.

г. Select an appropriate size PTA balloon catheter and dilate the lesion with conventional technique. The inflation diameter of the PTA balloon used for post dilatation should approximate the diameter of the reference vessel. Remove the PTA balloon from the patient.

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obtained.

b. Remove the guidewire and sheath and establish hemostasis. Note: Physician experience and discretion will determine the appropriate post-procedure drug regimen for each patient.

c. Close the entry wound as appropriate.

d. Discard the delivery system, guidewire, and sheath.

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The stent is MRI safe and does not interfere with, nor is affected by, the operation of an MRI device.

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• Appropriate antiplatelet/anticoagulant therapy should be administered pre and post procedure.

• The long term outcome following repeat dilatation of endothelialized stents is unknown at present.

• Persons allergic to nitinol (nickel titanium) may suffer an allergic reaction to this implant.

• Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

• The black dotted pattern on the grey temperature exposure indicator, found on the inner pouch, must be clearly visible. Do not use if entire temperature exposure indicator is completely black as the unconstrained stent diameter may have been compromised.

• Stenting across a major bile duct branch could lead to compromised future diagnostic or therapeutic procedures.

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- Any secondary expansion of the biliary stent should be performed with a device indicated for biliary stent placement.

- Do not use with Ethiodol or Lipiodol (Ethiodol and Lipiodol are trademarks of Guerbet S.A.) contrast media.

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- This product should only be used by physicians trained and experienced in diagnostic and interventional techniques. Standard techniques for interventional procedures should be employed.

- Do not use power injection systems with the delivery system.
- Use in patients with a history of contrast sensitivity is not recommended unless the patient can be adequately pre-medicated.
- Use of this device requires fluoroscopic control.
- The effect of heating in the MRI environment for overlapping stents or stents with fractured struts is not known.
- Safety and effectiveness has not been demonstrated in:
 - Lesions that are either totally or densely calcified
 - Pediatric patients

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System Handling - Precautions

- Store at ambient room conditions out of direct sunlight.
- Do not use beyond the "Use By" date.
- Carefully inspect the sterile packaging and device prior to use. Do not use if it appears damaged.
- Do not expose the delivery system to organic solvents.
- Do not attempt to remove the stent from the delivery system before use.

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- Do not use if the outer sheath of the device cannot be flushed.
- Do not use if the inner lumen of the guide-wire tube cannot be flushed.

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Stent Placement - Precautions

- If resistance is encountered at any time during the insertion procedure, do not force passage. Resistance may cause damage to the stent or lumen.
- If resistance occurs during movement through the sheath, carefully withdraw the stent system.
- If resistance is felt when initially retracting the outer deployment sheath, do not force deployment. Carefully withdraw the stent system without deploying the stent.
- Once stent deployment has begun, the stent must be fully deployed. The system is not designed for stent repositioning or recapturing.
- Use caution when crossing a deployed stent with any adjunct device.
- Prior to stent deployment, remove all slack from the catheter delivery system.

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There is no express or implied warranty, including without limitation any implied

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warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

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Cordis medical devices should be stored at temperatures ranging from +15°C to + 25°C, in a place protected from direct sunlight. Cordis medical devices being in operation are

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stable in temperature conditions from +32 °C to +42 °C and to biological liquids.

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Stents

PALMAZ GENESIS

The PALMAZ GENESIS Peripheral Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is supplied in four nominal lengths: 12 mm, 15 mm, 18 mm, 24 mm. Specific expansion ranges for each stent length are shown in Table 1. The PALMAZ GENESIS Peripheral Stent is sold unmounted for use with all catheter lengths of the Cordis OPTAR PRO Percutaneous Transluminal Angioplasty (PTA) Balloon Catheter [a .035" (.89 mm) diameter guidewire compatible system]. For further information see Table 1.

1. Use Cordis introducer tube INTR6 for all product codes and lengths.
2. Use of guiding catheters other than those recommended must have an inner lumen diameter of .088" (2.2 mm) or larger for 8F (2.7 mm) guides, .098" (2.5 mm) or larger for 9F (3.0 mm) guides, and .110" (2.8 mm) or larger for 10F (3.3 mm) guides.
3. These data are based on *in-vitro* testing.
4. Catheter suffix (*) refers to the usable catheter length, OPTA PRO is offered in the following lengths: T (40 cm), V (65 cm), S (80 cm), L (110 cm), X (135 cm). **Note:** Not all balloon sizes are sold in all catheter lengths.

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The PALMAZ GENESIS Peripheral Stent is indicated for use in the treatment of atherosclerotic disease of peripheral arteries below the aortic arch and for palliation of malignant neoplasms in the biliary tree.

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- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- A lesion that is within an aneurysm or an

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aneurysm with a proximal or distal segment adjacent to the lesion.

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- Embolization of atherosclerotic or thrombotic material
- Acute or sub-acute stent thrombosis
- Amputation
- Arteriovenous fistula
- Emergency surgery to correct vascular complications
- GI bleeding from anticoagulation/antiplatelet medication
- Hemorrhage/hematoma
- Artery injury, including perforation and dissection
- Rupture of retroperitoneum or neighboring organ
- Intimal tear

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- Pseudoaneurysm formation
- Restenosis of the stented artery
- Thrombus
- Tissue necrosis
- Total occlusion

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The patient may be started on enteric coated or nonenteric-coated aspirin 81-325 mg daily, one or two days prior to the procedure if deemed appropriate by the physician.

The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or sus-

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предыдущего, тромболизис должен precede стент deployment using standard acceptable practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

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4. Preparation of Stent Delivery Catheter - See Instructions for Use provided with the recommended balloon catheter.

a. Remove the balloon protector. b. Attach a stopcock to the catheter's inflation port. c. Open the stopcock and induce negative

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pressure.

d. Inflate the balloon by holding the proximal end of the catheter above the distal end. Hold the balloon vertically, with the balloon tip pointing down.

Warning: Inflation at a high rate may damage the balloon.

e. Continue pointing the distal tip of the balloon downward, and deflate the balloon again by inducing negative pressure. f. Ensure that the air in the balloon and the inflation lumen is removed. Repeat steps d–f. g. While holding negative pressure, close the stopcock. Maintain negative pressure until the stent is mounted. h. Looking proximal to distal, manually refold the balloon clockwise around the catheter.

5. Preparation of Stent Delivery System

a. Slide the stent over the distal end of the balloon, maintaining the balloon fold, until the radiopaque marker bands are equidistant from the ends of the stent.

b. Crimp the stent onto the balloon using the thumb and index finger of both hands. Apply even, concentric pressure along the length of the stent while fully rotating the assembly. Ensure that the stent is securely adhered to the balloon.

c. Open the stopcock allowing the inflation lumen and balloon to fill with diluted contrast medium.

d. Induce negative pressure to check the balloon for leaks, and then slowly release the negative pressure. e. Slightly pre-inflate both ends of balloon. Apply a pressure of no greater than one atmosphere.

f. Verify both balloon ends are slightly inflated.

Caution: Ensure that balloon inflation does not cause interference with guiding catheter or CSI.

g. Maintain low pressure during insertion of the stent/delivery assembly.

6. Stent Deployment Procedure for Peripheral Artery and Transhepatic Biliary Stent-

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ing

Insertion of Cordis Catheter Sheath Introducer (CSI) and Guidewire

а. Gain access at the appropriate site utilizing the recommended CSI size (see Table 1).

Caution: Always use a CSI for the implant procedure to protect the puncture site and, in the case of biliary stenting, the liver tract, and to avoid dislodging the stent from the balloon.

б. Insert the appropriately sized guidewire that is at least twice the length of the delivery system across the lesion to be stented through the CSI.

Note: If the physician determines that predilation is necessary, PTA techniques may be used. Maintain lesion access with the guidewire.

Caution: Avoid overexpansion of the balloon to prevent bleeding, dissection or patient discomfort.

с. Advance the CSI completely across the lesion.

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а. Flush guidewire lumen of the delivery system. Position the flared end of the introducer tube over the distal tip of the stented balloon and slide forward to protect the stent during CSI insertion. Backload onto the guidewire. Place the assembly through the CSI hemostatic valve until resistance is met. Carefully advance the stent and delivery system through the introducer tube and hemostatic valve. When the stent has passed into the body of the CSI, remove the introducer tube.

б. Continue to advance the device over the guidewire through the hemostatic valve and CSI.

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Stent Deployment

а. Keeping the CSI immobile, observe fluoroscopically as the stent is advanced through the CSI to the site of the lesion.

Caution: If excessive resistance is encountered, do not force passage.

б. Under fluoroscopy, use the balloon marker bands and the radiopaque stent to position the stent centrally within the lesion. During positioning, verify that the stent is still centered within the balloon marker bands and has not been dislodged.

с. Prior to stent expansion, ensure that stent and balloon are completely exposed from the CSI or guiding catheter.

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д. Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the catheter label. Expand the diameter of the stent to the diameter of the reference vessel.

Caution: Do not exceed the rated burst pressure shown on the catheter label.

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Magnetic Resonance Imaging (MRI) Compatibility

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7. Delivery System Withdrawal

- a. After deploying the stent, deflate the balloon by pulling a vacuum, allowing adequate time for the balloon to fully deflate prior to removal.
- b. Carefully rotate balloon counterclockwise to ensure separation of the balloon from the stent.
- c. While maintaining negative pressure on the balloon, slowly withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.
- d. Withdraw the deflated delivery system into the CSI or guiding catheter.
- e. Remove and discard the delivery system, guidewire and CSI or guiding catheter from the body.

The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the stainless steel stent should not be greater than those caused by metal surgical clips. In a study of larger size stainless steel stents, no stent migration was observed in an MRI field strength up to 1.5 Tesla. However, an MRI scan should not be performed until the stent implantation site has had a chance to heal (estimated to be 8 weeks), in order to further minimize the risk of migration.

- Do not use if the inner package is opened or damaged.
- The **PALMAZ GENESIS** Peripheral Stent is intended for single use only. Do not resterilize and/or reuse the device.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and resterilization after initial use.
- Reuse of this product, including after reprocessing and/or re-sterilization, may

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cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

- Store in a cool, dark, dry place.
- Use the stent prior to the "Use By" date specified on the package.
- Persons with allergic reactions to stainless steel or its components (for example, nickel) may suffer an allergic response to this implant.
- Once fully deployed, the stent cannot be repositioned.
- Stenting across a bifurcation could compromise future diagnostic or therapeutic procedures.
- Avoid stent placement which would obstruct access to a vital side branch.
- It is recommended that stents not be implanted in patients with confirmed pregnancies.

For Peripheral Artery Stent Implantation Only:

- The safety and efficacy of the **PALMAZ GENESIS** Peripheral Stent have not been established for use in patients for the treatment of atherosclerotic disease of the femoral artery and other superficial arteries.
- As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or the retroperitoneum.
- The stent may cause a thrombus, distal embolization or may migrate from the site of the implant down the arterial lumen.
- Overstretching of the artery may result in rupture and life threatening bleeding.
- When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.

The device should only be used by physicians who are trained in interventional techniques such as percutaneous

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transluminal angioplasty, placement of stents, and transhepatic access.

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- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception, stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.

- To assure full expansion, inflate to at least the nominal pressure recommended on the catheter label.

- Do not exceed the rated burst pressure recommended on the catheter label.

- Ensure stent is securely crimped onto balloon prior to use.

- Do not attempt to remove or readjust the stent once crimped onto the delivery system.

- When catheters are in the body, they should be manipulated only under fluoroscopy.

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- When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chance for dislodging the proximal stent.

- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.

- In the event of complications such as infections, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.

- Caution should be taken when stenting patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.

- In patients requiring the use of antacids and/or H₂-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g. aspirin) may be adversely affected.

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Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential

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damages resulting from reuse of the product.

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Instruction for Use Stents

PALMAZ GENESIS

The PALMAZ GENESIS Peripheral Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is supplied in five nominal lengths: 19 mm, 25 mm, 29 mm, 39 mm and 59 mm. Specific expansion ranges for each stent length are shown in Table 1. The PALMAZ GENESIS Peripheral Stent is sold unmounted for use with all catheter lengths of the Cordis OPTA PRO Percutaneous Transluminal Angioplasty (PTA) Balloon Catheter (a .035" [0.89 mm] diameter guidewire compatible system). For further information, see Table 1.

1. Use Cordis introducer tube INTR6 for all product codes and lengths.

2 Use of guiding catheters other than those recommended must have an inner lumen diameter of .110" (2.8 mm) or larger for 10F (3.3 mm) guides.

3 These data are based on *in-vitro* testing.

4 Catheter suffix (*) refers to the usable catheter length. OPTA PRO is offered in the following lengths: T (40 cm), V (65 cm), S (80 cm), L (110 cm), X (135 cm).

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Generally, contraindications to percutaneous transluminal angioplasty (PTA) are al-

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so contraindications for stent placement. Contraindications include, but may not be limited to:

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Potential complications associated with **peripheral artery stent** implantation may include, but are not limited to, the follow-

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ing:

- Embolization of atherosclerotic or thrombotic material
- Acute or sub-acute stent thrombosis
- Amputation
- Arteriovenous fistula
- Emergency surgery to correct vascular complications
- GI bleeding from anticoagulation/antiplatelet medication
- Hemorrhage/hematoma
- Artery injury, including perforation and dissection
- Rupture of retroperitoneum or neighboring organ
- Intimal tear
- Pseudoaneurysm formation
- Restenosis of the stented artery
- Thrombus
- Tissue necrosis
- Total occlusion

Potential complications associated with **transhepatic biliary endoprostheses** implantation may include, but are not limited to, the following:

- Stent obstruction secondary to tumor growth through the stent
- Tumor overgrowth at the stent ends
- Bile duct occlusion/obstruction
- Pancreatitis
- Bile duct perforation
- Parenchymal hemorrhage
- Hemobilia
- Peritonitis
- Abscess
- Rupture, overstretching of duct
- Cholangitis

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2. Select Stent Size

а. Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent (see Table 1).

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б. Measure the diameter of the reference vessel to determine the appropriate size stent and delivery system.

3. Preparation of Stent

а. Open the carton to reveal the pouch containing the stent.

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lumen and balloon to fill with diluted contrast medium.

д. Induce negative pressure to check the balloon for leaks, and then slowly release the negative pressure. е. Slightly pre-inflate both ends of balloon. Apply a pressure of no greater than one atmosphere.

ф. Verify both balloon ends are slightly inflated.

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г. Maintain low pressure during insertion of the stent/delivery assembly.

6. Stent Deployment Procedure for Peripheral Artery and Transhepatic Biliary Stenting

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The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the stainless steel stent should not be greater than those caused by metal surgical clips. In a study of larger size stainless steel stents, no stent migration was observed in an MRI field strength up to 1.5 Tesla. However, an MRI scan should not be performed until the stent implantation site has had a chance to heal (estimated to be 8 weeks), in order to further mini-

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- The safety and efficacy of the PALMAZ GENESIS Peripheral Stent have not been established for use in patients for the treatment of atherosclerotic disease of the femoral artery and other superficial arteries.
- As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or the retroperitoneum.
- The stent may cause a thrombus, distal embolization or may migrate from the site of the implant down the arterial lumen.
- Overstretching of the artery may result in rupture and life threatening bleeding.
- When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.

- The device should only be used by physicians who are trained in interventional techniques such as percutaneous transluminal angioplasty, placement of stents, and transhepatic access.
- The delivery system is not designed for use with power injection systems.
- Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.
- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception, stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.
- To assure full expansion, inflate to at least the nominal pressure recommended on the

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catheter label.

- Do not exceed the rated burst pressure recommended on the catheter label.
- Ensure stent is securely crimped onto balloon prior to use.
- Do not attempt to remove or readjust the stent once crimped onto the delivery system.
- When catheters are in the body, they should be manipulated only under fluoroscopy.

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- When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chance for dislodging the proximal stent.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infections, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.
- Caution should be taken when stenting patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.
- In patients requiring the use of antacids and/or H₂-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g. aspirin) may be adversely affected.

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Authorized Representative of the manufacturer in the Russian Federation:

Johnson & Johnson LLC

Russian Federation, Krylatskaya street 17, building 2, Moscow, 121614

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Instruction for Use Stents PALMAZ GENESIS

The PALMAZ GENESIS Peripheral Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is supplied in four nominal lengths: 29 mm, 39 mm, 59 mm and 79 mm. Specific expansion ranges for each stent length are shown in Table 1. The PALMAZ GENESIS Peripheral Stent is sold unmounted for use with all catheter lengths of the Cordis OPTA PRO Percutaneous Transluminal Angioplasty (PTA) Balloon Catheter [a .035" (.89 mm) diameter guidewire compatible system]. For further information, see Table 1.

1. Use Cordis introducer tube INTR6 for all product codes and lengths.
2. Use of guiding catheters other than those recommended must have an inner lumen diameter of .088" (2.2 mm) or larger for 8F (2.7 mm) guides, .098" (2.5 mm) or larger for 9F (3.0 mm) guides, and .110" (2.8 mm) or larger for 10F (3.3 mm) guides.
3. These data are based on in-vitro testing.
4. Catheter suffix (*) refers to the usable catheter length. OPTA PRO is offered in the following lengths: T (40 cm), V (65 cm), S (80 cm), L (110 cm), X (135 cm).

Note: Not all balloon sizes are sold in all catheter lengths.

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- Patients with highly calcified lesions resistant to PTA.

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- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- A lesion that is within an aneurysm or an aneurysm with a proximal or distal segment adjacent to the lesion.

Transhepatic Biliary Stent Implantation

Contraindications associated with the use of transhepatic biliary endoprotheses include:

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Potential complications associated with **peripheral artery and transhepatic biliary stent** implantation may include, but are not limited to, the following:

- Sepsis/infection
- Stent migration/embolization
- Drug reaction, allergic reaction to contrast medium.

Potential complications associated with **peripheral artery stent** implantation may include, but are not limited to, the following:

- Embolization of atherosclerotic or thrombotic material
- Acute or sub-acute stent thrombosis
- Amputation
- Arteriovenous fistula

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- Emergency surgery to correct vascular complications
- GI bleeding from anticoagulation/antiplatelet medication
- Hemorrhage/hematoma
- Artery injury, including perforation and dissection
- Rupture of retroperitoneum or neighboring organ
- Intimal tear
- Pseudoaneurysm formation
- Restenosis of the stented artery
- Thrombus
- Tissue necrosis
- Total occlusion

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The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

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2. Select Stent Size

a. Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent (see Table 1).

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b. Measure the diameter of the reference vessel to determine the appropriate size stent and delivery system.

3. Preparation of Stent

a. Open the carton to reveal the pouch containing the stent.

b. Inspect the stent package for damage to the sterile barrier. Remove the stent from the package and rinse in sterile heparinized saline.

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4. *Preparation of Stent Delivery Catheter - See Instructions for Use provided with the recommended balloon catheter.*

a. Remove the balloon protector. b. Attach a stopcock to the catheter's inflation port. c. Open the stopcock and induce negative pressure.

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e. Continue pointing the distal tip of the balloon downward, and deflate the balloon again by inducing negative pressure. f. Ensure that the air in the balloon and the inflation lumen is removed. Repeat steps d-f. g. While holding negative pressure, close the stopcock. Maintain negative pressure until the stent is mounted. h. Looking proximal to distal, manually refold the balloon clockwise around the catheter.

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d. Induce negative pressure to check the balloon for leaks, and then slowly release the negative pressure. e. Slightly pre-inflate both ends of balloon. Apply a pressure of

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the guidewire. Place the assembly through the CSI hemostatic valve until resistance is met. Carefully advance the stent and delivery system through the introducer tube and hemostatic valve. When the stent has passed into the body of the CSI, remove the introducer tube.

b. Continue to advance the device over the guidewire through the hemostatic valve and CSI.

Stent Deployment

a. Keeping the CSI immobile, observe fluoroscopically as the stent is advanced through the CSI to the site of the lesion.

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b. Under fluoroscopy, use the balloon marker bands and the radiopaque stent to position the stent centrally within the lesion. During positioning, verify that the stent is still centered within the balloon marker bands and has not been dislodged.

c. Prior to stent expansion, ensure that stent and balloon are completely exposed from the CSI or guiding catheter.

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d. Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the catheter label. Expand the diameter of the stent to the diameter of the reference vessel.

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e. Remove and discard the delivery system, guidewire and CSI or guiding catheter from the body.

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- The device should only be used by physicians who are trained in interventional techniques such as percutaneous transluminal angioplasty, placement of stents, and transhepatic access.
- The delivery system is not designed for use with power injection systems.
- Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.
- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception, stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.
- To assure full expansion, inflate to at least the nominal pressure recommended on the catheter label.

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Instruction for Use Stents PALMAZ GENESIS on SLALOM

The medical device is the PALMAZ GENESIS Peripheral Stent on SLALOM .018" Delivery System. For stent and balloon dimensions, see the carton label.

Stent

The PALMAZ GENESIS Peripheral Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is sold mounted on a Cordis percutaneous transluminal angioplasty (PTA) balloon catheter.

Delivery System

A Cordis PTA balloon catheter is used as a delivery system for the PALMAZ GENESIS Peripheral Stent. The delivery system utilizes an over-the-wire design and is a catheter with a distal DURALYN balloon and a proximal Y-connector. Two radiopaque marker bands aid in stent placement. The injectate lumen (marked "THRU") is used to inject contrast medium via hand injection, and to track the catheter over a guidewire. The inflation lumen (marked "BALLOON") is used to inflate and deflate the balloon. The nominal balloon size is printed on the Y-connector.

Introducer Tube

A metal introducer tube is included in the packaging.

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The Cordis PALMAZ GENESIS Peripheral Stent on SLALOM .018" Delivery System is intended for use in the treatment of atherosclerotic disease of peripheral arteries below the aortic arch and for palliation of malignant neoplasms in the biliary tree

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- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- A lesion that is within an aneurysm or an aneurysm with a proximal or distal segment adjacent to the lesion.

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- Stenting of a perforated duct where leakage from the duct could be exacerbated by the prosthesis.
- Patients with uncorrected bleeding disorders.
- Severe ascites.

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Potential complications associated with **peripheral artery and transhepatic biliary stent implantation** may include, but are not limited to, the following:

- Sepsis/infection
- Stent migration.

Potential complications associated with **peripheral artery stent implantation** may include, but are not limited to, the following:

- Embolization of atherosclerotic or thrombotic material
- Acute or sub-acute stent thrombosis
- Amputation
- Arteriovenous fistula
- Emergency surgery to correct vascular complications
- GI bleeding from anticoagulation/antiplatelet medication
- Hemorrhage/hematoma
- Artery injury, including perforation and dissection
- Rupture of retroperitoneum or neighboring organ
- Intimal tear
- Pseudoaneurysm formation
- Restenosis of the stented artery
- Thrombus
- Tissue necrosis
- Total occlusion
- Stent Embolisation
- Drug reaction, allergic reaction to contrast medium
- Air embolism.

Potential complications associated with **transhepatic biliary endoprostheses implantation** may include, but are not limited to, the following:

- Stent obstruction secondary to tumor growth through the stent
- Tumor overgrowth at the stent ends
- Bile duct occlusion/obstruction
- Pancreatitis
- Bile duct perforation
- Parenchymal hemorrhage
- Hemobilia
- Peritonitis
- Abscess

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- Rupture, overstretching of duct
- Cholangitis.

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1. The patient may be started on enteric coated or nonenteric-coated aspirin 81-325 mg daily, one or two days prior to the procedure if deemed appropriate by the physician.
2. The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

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1. Perform standard diagnostic angiography or percutaneous cholangiogram, as applicable, to evaluate the lesion.

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2. Selection of Stent Size

a. Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.

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g. Hold the syringe and proximal end of the catheter above the distal end of the catheter, and hold the balloon vertically with the balloon tip pointing down.

h. Close the stopcock to the inflation port. Remove the syringe.

i. To ensure air contained in the balloon and inflation lumen is removed, apply negative pressure twice as instructed in steps f-h.

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Caution: Never re-advance a CSI or guiding catheter over an exposed stent to avoid dislodging the stent.

Caution: If a Tuohy-Borst type adjustable hemostasis valve is used, avoid over tightening since

this may restrict the flow of contrast media in and out of the balloon, thereby slowing inflation/deflation.

d. Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the label. Expand the diameter of the stent to the diameter of the reference vessel.

Caution: Do not exceed the rated burst pressure shown on the label.

Note: It is strongly recommended that the guidewire, the delivery system, or both, remain across the lesion until the procedure is completed and the delivery system is removed from the vessel.

Note: For peripheral artery stenting, if the patient has significant femoral artery disease or if additional surgical intervention is contemplated, the stent may be placed in an operating room following surgical cut-down and isolation of the femoral artery.

5. Delivery System Withdrawal

a. After deploying the stent, deflate the balloon by pulling a vacuum, allowing adequate time for the balloon to fully deflate prior to removal.

b. Carefully rotate balloon counterclockwise to ensure separation of the balloon from the stent.

c. While maintaining negative pressure on the balloon, slowly withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.

d. Withdraw the deflated delivery system into the CSI or guiding catheter.

e. Remove and discard the delivery system, guidewire and CSI and if applicable guiding catheter from the body.

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The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the stainless steel stent should not be greater than those caused by metal surgical clips. In a study of larger size stainless steel stents, no stent migration was observed in an MRI field strength up to 1.5 Tesla.* However, an MRI scan should not be performed until the stent implantation site has had a chance to heal (estimated to be 8 weeks), in order to further minimize the risk of migration.

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Never use air or any gaseous medium to inflate the balloon.

- Do not use Ethiodol** (Ethiodol is a trademark of Guerbet S.A.) or Lipiodol contrast medium.
- Do not expose the delivery system to organic solvents (e.g. alcohol).
- Do not autoclave. Exposure to temperatures above 54°C (130°F) may damage the catheter.
- Persons with allergic reactions to stainless steel may suffer an allergic response to this implant.
- Once fully deployed, the stent cannot be repositioned.
- Stenting across a bifurcation could compromise future diagnostic or therapeutic procedures.

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- The safety and efficacy of the PALMAZ GENESIS Peripheral Stent on SLALOM .018" Delivery System have not been established for use in patients for treatment of atherosclerotic disease of the femoral artery.
- As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or the retroperitoneum.
- The stent may cause a thrombus, distal embolization or may migrate from the site of the implant down the arterial lumen.
- Overstretching of the artery may result in rupture and life threatening bleeding.
- When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.
- Avoid stent placement which would ob-

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struct access to a vital side branch.

- It is recommended that stents not be implanted in patients with confirmed pregnancies.
- To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

- The device should only be used by physicians who are trained in interventional techniques such as PTA, placement of stents and transhepatic access.
- Prior to stenting, the PALMAZ GENESIS Peripheral Stent on SLALOM .018" Delivery System should be examined to verify functionality and integrity.
- The delivery system is not designed for use with power injection systems.
- Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.
- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception: stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.
- Do not attempt to remove or readjust the stent on the delivery system.
- To assure full expansion, inflate to at least the recommended nominal pressure as shown on the label.
- Do not exceed the rated burst pressure recommended on the label. The rated burst pressure is based on the results of in vitro testing. At least 99.9% of the balloons (with a 95% confidence), will not burst at or below

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their rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization.

- When catheters are in the body, they should be manipulated only under fluoroscopy. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum.
- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the delivery system through a smaller size sheath introducer than indicated on the label.

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- When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chance for dislodging the proximal stent.
- In the event of thrombosis of the expanded stent, thrombolysis and PTA should be attempted.
- In the event of complications such as infections, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.
- Caution should be taken when stenting patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.
- In patients requiring the use of antacids and/or H2-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g. aspirin) may be ad-

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versely affected.

- Not intended for precise arterial blood pressure monitoring.
- Consider the use of systemic heparinization. Flush all devices entering vascular system with sterile heparinized saline or similar isotonic solution.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do

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not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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Instruction for Use Stents PALMAZ GENESIS ON OPTA PRO

For stent and balloon dimensions, see the carton label.

Stent

The PALMAZ GENESIS Transhepatic Biliary Stent is a balloon-expandable, laser cut stent made from 316L stainless steel tubing. The stent is sold mounted on a Cordis balloon catheter.

Delivery System

A Cordis balloon catheter is used as a delivery system for the PALMAZ GENESIS Transhepatic Biliary

Stent. The delivery system utilizes an over-the-wire design and is a catheter with a distal DURALYN

balloon and a proximal Y-connector. Two radiopaque marker bands aid in stent placement. The injectate lumen (marked "THRU") is used to inject contrast medium via hand injection, and to track the catheter over a guidewire. The inflation lumen (marked "BALLOON") is used to inflate and deflate the balloon. The nominal balloon size is printed on the Y-connector.

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The Cordis PALMAZ GENESIS Peripheral Stent on OPTA PRO .035" Delivery System is not indicated for use in the arcus aorta or the aorta descendens to the bifurcation aortae.

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Peripheral Artery Stent Implantation

Generally, contraindications to percutaneous transluminal angioplasty (PTA) are also contraindications for stent placement. Contraindications include, but may not be limited to:

Patients with highly calcified lesions resistant to PTA.

Patients with a target lesion with a large amount of adjacent acute or subacute thrombus.

Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.

Patients with perforated vessels evidenced by extravasation of contrast media.

A lesion that is within an aneurysm or an aneurysm with a proximal or distal segment adjacent to the lesion.

Transhepatic Biliary Stent Implantation

Contraindications associated with the use of transhepatic biliary endoprostheses include :

Stenting of a perforated duct where leakage from the duct could be exacerbated by the prosthesis.

Patients with uncorrected bleeding disorders. Severe ascites.

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Potential complications associated with peripheral artery and transhepatic biliary stent implantation may include, but are not limited to, the following:

Sepsis/infection
Stent migration.

Potential complications associated with peripheral artery stent implantation may include, but are not limited to, the following:

- Aneurysm
- Allergic reaction to stent material cobalt chromium or one of its components, or to contrast medium
- Arteriovenous fistula
- Cerebrovascular accident/Stroke
- Death
- Embolization (air, tissue, stent, atherosclerotic or thrombotic material, cholesterol)
- Emergency surgery to correct vascular complications (revascularization/bypass)
- Failure to deliver the stent to the intended site
- Hemorrhage/hematoma
- Hypertension, Hypotension
- Infection/inflammation
- Myocardial ischemia/infarction and arrhythmia
- Peripheral neuropathy
- Persistent abdominal pain
- Renal insufficiency, renal failure/dialysis, renal infarction, renal transplantation/Nephrectomy
- Restenosis of the stented artery
- Rupture of retroperitoneum or neighboring organ
- Tissue necrosis/ulceration
- Stent migration/embolization
- Stent thrombosis/occlusion
- Vascular complications (e.g. at puncture and/or lesion site, dissection, intimal tear, pseudoaneurysm, rupture and perforation, spasm, occlusion).

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Potential complications associated with transhepatic biliary endoprosthesis implantation may include, but are not limited to, the following:

- Stent obstruction secondary to tumor growth through the stent
- Tumor overgrowth at the stent ends
- Bile duct occlusion/obstruction
- Pancreatitis
- Bile duct perforation
- Parenchymal hemorrhage
- Hemobilia
- Peritonitis
- Abscess
- Rupture, overstretching of duct
- Cholangitis.

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Pre-Procedure for Peripheral Artery Stenting

- The patient may be started on enteric coated or nonenteric-coated aspirin 81-325 mg daily, one or two days prior to the procedure if deemed appropriate by the physician.
- The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed

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to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

Procedure for Peripheral Artery and Transhepatic Biliary Stenting

- Perform standard diagnostic angiography or percutaneous cholangiogram, as applicable, to evaluate the lesion.
- *Selection of Stent Size*
Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.

Note: Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.

- Measure the diameter of the reference vessel to determine the appropriate size stent and delivery system

Preparation of Stent Delivery System

1. Open the outer box and reveal the inner package containing the tray with the stent and delivery system and introducer tube.

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Deployment Procedure for Peripheral Artery and Transhepatic Biliary Stenting

- Insertion of Cordis Sheath Introducer (CSI) and Guidewire
- Gain access at the appropriate site utilizing the CSI size recommended on the label. **Caution:** Always use a CSI for the implant procedure to protect the

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puncture site and, in the case of biliary stenting, the liver tract, and to avoid dislodging the stent from the balloon.

Insert a guidewire that is at least twice the length of the delivery system across the lesion to be stented through the CSI. **Note:** Compatible guidewire diameter is listed on the label. **Note:** If the physician determines that predilation is necessary, standard PTA techniques may be used. Maintain lesion access with the guidewire. **Caution:** Avoid overexpansion of the balloon to prevent bleeding, dissection or patient discomfort.

For transhepatic biliary stent implantation only:

Advance the CSI completely across the lesion.

Optional: Select and insert the appropriate size guiding catheter if the CSI utilized during the procedure cannot be advanced to the lesion site.

Introduction of Stent Delivery System

Flush guidewire lumen of the delivery system. Position the flared end of the introducer tube over the distal tip of the stented balloon and slide forward to protect the stent during CSI insertion. Backload onto the guidewire. Place the assembly through the CSI hemostatic valve until resistance is met. Carefully advance the stent and delivery system through the introducer tube and hemostatic valve. When the stent has passed into the body of the CSI, remove the introducer tube. Continue to advance the device over the guidewire through the hemostatic valve and CSI.

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Stent Deployment

Keeping the CSI immobile, observe fluoroscopically as the stent is advanced through the CSI to the site of the lesion. **Caution:** If strong resistance is met during advancement or withdrawal of the catheter, discontinue movement and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system.

Under fluoroscopy, use the balloon marker bands and the radiopaque stent to position the stent centrally within the lesion. During positioning, verify that the stent is still centered within the balloon marker bands and has not been dislodged.

Prior to stent expansion, ensure that stent and balloon are completely exposed from the CSI or guiding catheter. **Caution:** Never re-advance a CSI or guiding catheter over an exposed stent to avoid dislodging the stent. **Caution:** If a Tuohy-Borst type adjustable hemostasis valve is used, avoid over tightening since this may restrict the flow of contrast media in and out of the balloon, thereby slowing inflation/deflation.

Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the label. Expand the diameter of the stent to the diameter of the reference vessel. **Caution:** Do not exceed the rated burst pressure shown on the label. **Note:** It is strongly recommended that the guidewire, the delivery system, or both, remain across the lesion until the procedure is completed and the delivery system is removed from the vessel. **Note:** For peripheral artery stenting, if the patient has significant femoral artery disease or if additional surgical intervention is contemplated, the stent may be placed in an operating room following surgical cut-down and isolation of the femoral ar-

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Delivery System Withdrawal

- After deploying the stent, deflate the balloon by pulling a vacuum, allowing adequate time for the balloon to fully deflate prior to removal.
- Carefully rotate balloon counterclockwise to ensure separation of the balloon from the stent.
- While maintaining negative pressure on the balloon, slowly withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.
- Withdraw the deflated delivery system into the CSI or guiding catheter.
- Remove and discard the delivery system, guidewire and CSI and if applicable guiding catheter from the body.

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For Peripheral Artery Stent Implantation Only

The safety and efficacy of the PALMAZ GENESIS Peripheral Stent on OPTA PRO .035" Delivery System have not been estab-

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lished for use in patients for treatment of atherosclerotic disease of the femoral artery.

As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or the retroperitoneum.

The stent may cause a thrombus, distal embolization or may migrate from the site of the implant down the arterial lumen.

Overstretching of the artery may result in rupture and life threatening bleeding.

When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.

Avoid stent placement which would obstruct access to a vital side branch.

It is recommended that stents not be implanted in patients with confirmed pregnancies.

To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis

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The device should only be used by physicians who are trained in interventional techniques such as PTA, placement of stents and transhepatic access.

- Prior to stenting, the PALMAZ GENESIS Peripheral Stent on OPTA PRO .035" Delivery System should be examined to verify functionality and integrity.

- The delivery system is not designed for use with power injection systems.

- Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.

- To avoid the possibility of dissimilar metal

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corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception: stents made of 316L stainless steel are compatible with stents made of nickel titanium alloy.

- Do not attempt to remove or readjust the stent on the delivery system.

To assure full expansion, inflate to at least the recommended nominal pressure as shown on the label.

- Do not exceed the rated burst pressure recommended on the label. The rated burst pressure is based on the results of *in-vitro* testing. At least 99.9% of the balloons (with a 95% confidence), will not burst at or below their rated burst pressure. Use of a pressure monitoring device is recommended to prevent over pressurization.
- When catheters are in the body, they should be manipulated only under fluoroscopy. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum.
- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the delivery system through a smaller size sheath introducer than indicated on the label.

For Peripheral Artery Stent Implantation Only

When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent in placement of the distal stent and reduces the chance for dislodging the proximal stent.

In the event of thrombosis of the expanded

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stent, thrombolysis and PTA should be attempted.

In the event of complications such as infections, pseudoaneurysm or fistulization, surgical removal of the stent may be required. Standard surgical procedure is appropriate.

Caution should be taken when stenting patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.

In patients requiring the use of antacids and/or H₂-antagonists before or immediately after stent placement, oral absorption of antiplatelet agents (e.g. aspirin) may be adversely affected. Not intended for precise arterial blood pressure monitoring.

Consider the use of systemic heparinization. Flush all devices entering vascular system with sterile heparinized saline or similar isotonic solution

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The manufacturer warrants conformance of this medical device with the declared characteristics within shelf life (sterility term) provided that all requirements of the manufacturer's use instructions are met.

There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty

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except as specifically set forth herein.

Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for any direct, incidental, or consequential damages resulting from reuse of the product.

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**Instruction for Use
Stents
PALMAZ BLUE**

The medical device is Peripheral Stent PALMAZ BLUE on SLALOM .018" Delivery System. For stent and balloon dimensions, see the carton label.

Stent

The PALMAZ BLUE Peripheral Stent is a balloon-expandable, laser cut stent made from cobalt-chromium tubing. The stent is sold mounted on a Cordis percutaneous transluminal angioplasty (PTA) balloon catheter.

Delivery System

A Cordis PTA balloon catheter is used as a delivery system for the PALMAZ BLUE Peripheral Stent. The delivery system utilizes an over-the-wire design and is a catheter with a distal DURALYN balloon and a proximal Y-connector. Two radiopaque marker bands aid in stent placement. The injectate lumen (marked "THRU") is used to inject contrast medium via hand injection, and to track the catheter over a guidewire. The inflation lumen (marked "BALLOON") is used to inflate and deflate the balloon. The nominal balloon size is printed on the Y-connector.

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The Cordis PALMAZ BLUE .018 Peripheral Stent System is intended for use in the treatment of atherosclerotic disease of peripheral arteries below the aortic arch and for palliation of malignant neoplasms in the biliary tree.

The Cordis PALMAZ BLUE .018 Peripheral Stent System is not indicated for use in the arcus aorta or the aorta descendens to the bifurcation aortae.

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Transhepatic Biliary Stent Implantation

Contraindications associated with the use of transhepatic biliary endoprostheses include:

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Patients with uncorrected bleeding disorders.

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Severe ascites.

Potential complications associated with **pe-ripheral artery and transhepatic biliary stent implantation** may include, but are not limited to, the following:

Sepsis/infection
Stent migration.

Potential complications associated with **pe-ripheral artery stent implantation** may in-clude, but are not limited to, the following:

- Aneurysm
- Allergic reaction to stent material cobalt chromium or one of its components, or to contrast medium
- Arteriovenous fistula
- Cerebrovascular accident/Stroke
- Death
- Embolization (air, tissue, stent, athero-sclerotic or thrombotic material, cholest-erol)
- Emergency surgery to correct vascular complications (revascularization/bypass)
- Failure to deliver the stent to the intend-ed site
- Hemorrhage/hematoma
- Hypertension, Hypotension
- Infection/inflammation
- Myocardial ischemia/infarction and ar-rhythmia
- Peripheral neuropathy
- Persistent abdominal pain
- Renal insufficiency, renal fail-ure/dialysis, renal infarction, renal trans-plantation/ Nephrectomy
- Restenosis of the stented artery
- Rupture of retroperitoneum or neighbor-ing organ
- Tissue necrosis/ulceration
- Stent migration/embolization
- Stent thrombosis/occlusion
- Vascular complications (e.g.at puncture and/or lesion site, dissection, intimal tear, pseudoaneurysm, rupture and perfo-ration, spasm, occlusion).

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1. The patient may be started on enteric coated or nonenteric-coated aspirin 81-325 mg daily, one or two days prior to the procedure if deemed appropriate by the physician.
2. The percutaneous placement of the stent in a stenotic or obstructed artery should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion(s) and the collateral flow. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard acceptable

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practice. Access vessels must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

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2. Selection of Stent Size

a. Measure the length of the target lesion or stricture to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion or stricture. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent.

Note: Should more than one stent be required, place the stent most distal from the puncture or incision site first, followed by placement of the proximal stent in tandem.

b. Measure the diameter of the reference vessel or duct to determine the appropriate size stent and delivery system.

3. Preparation of Stent Delivery System

a. Open the outer box and reveal the inner package containing the tray with the stent and delivery system and introducer tube.

b. Open the inner package and carefully extract the tray with the stent and delivery system and introducer tube.

c. Hold the tray in one hand. With the other hand, gently grasp the hub and carefully remove the stent/delivery system from the tray. Remove the introducer tube from the tray and discard the tray.

d. Inspect the crimped stent for adherence to the balloon and centered placement in relation to the balloon marker bands. Do not reposition the stent or hand crimp.

e. Attach a stopcock to the catheter's inflation port.

f. Attach a syringe, open the stopcock and induce negative pressure.

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catheter above the distal end of the catheter, and hold the balloon vertically with the balloon tip pointing down.

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j. Connect an angioplasty inflation system. Open the stopcock and slowly fill the inflation lumen and the balloon with diluted contrast medium.

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1. Insertion of Cordis Sheath Introducer (CSI) and Guidewire

a. Gain access at the appropriate site utilizing the CSI size recommended on the label.

Caution: Always use a CSI for the implant procedure to protect the puncture or incision site and, in the case of biliary stenting, the liver tract, and to avoid dislodging the stent from the balloon.

b. Insert a guidewire that is at least twice the length of the delivery system across the lesion or stricture to be stented through the CSI.

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Note: If the physician determines that predilation is necessary, standard dilatation techniques may be used. Maintain lesion or stricture access with the guidewire.

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tion only: Advance the CSI completely across the stricture.

2. Optional: Select and insert the appropriate size guiding catheter if the CSI utilized during the procedure cannot be advanced to the lesion site.

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a. Flush guidewire lumen of the delivery system. Position the flared end of the introducer tube over the distal tip of the stented balloon and slide forward to protect the stent during CSI insertion. Backload onto the guidewire. Place the assembly through the CSI hemostatic valve until resistance is met. Carefully advance the stent and delivery system through the introducer tube and hemostatic valve. When the stent has passed into the body of the CSI, remove the introducer tube.

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4. Stent Deployment

a. Keeping the CSI immobile, observe fluoroscopically as the stent is advanced through the CSI to the site of the lesion or stricture.

Caution: If strong resistance is met during advancement or withdrawal of the catheter, discontinue movement and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, withdraw the entire system.

b. Under fluoroscopy, use the balloon marker bands and the radiopaque stent to position the stent centrally within the lesion or stricture. During positioning, verify that the stent is still centered within the balloon marker bands and has not been dislodged.

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Caution: If a Tuohy-Borst type adjustable hemostasis valve is used, avoid over tightening since this may restrict the flow of contrast media in and out of the balloon,

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thereby slowing inflation/deflation.

- d. Using an inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the label. Expand the diameter of the stent to the diameter of the reference vessel or duct.

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Note: It is strongly recommended that the guidewire, the delivery system, or both, remain across the lesion or stricture until the procedure is completed and the delivery system is removed from the vessel or duct.

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5. Delivery System Withdrawal

- a. After deploying the stent, deflate the balloon by pulling a vacuum, allowing adequate time for the balloon to fully deflate prior to removal.

- b. Carefully rotate balloon counterclockwise to ensure separation of the balloon from the stent.

- c. While maintaining negative pressure on the balloon, slowly withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.

- d. Withdraw the deflated delivery system into the CSI or guiding catheter.

- e. Remove and discard the delivery system, guidewire and CSI and if applicable guiding catheter from the body.

Note: Packaging disposal can be handled per normal industrial standard. The used device shall be handled according hospital procedures

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The Cordis PALMAZ BLUE Stent has been shown to be MR safe, when tested either individually or in the combination of 2-stents overlapped with 4 mm, at field strengths of 3.0 Tesla or less, a maximum spatial gradient of 5 T/m, and a maximum whole body averaged specification absorption rate (SAR) of 4.0 W/kg for 15:00 min of MR imaging. MR imaging quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the PALMAZ BLUE Stent.

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tion medium (a 50/50 mixture by volume of contrast medium and normal saline). Never use air or any gaseous medium to inflate the balloon.

- Do not use Ethiodol or Lipiodol* contrast medium.

- Do not expose the delivery system to organic solvents (e.g. alcohol).

- Exposure to temperatures above 54°C (130°F) may damage the catheter.

- Persons with allergic reactions to a cobalt-chromium alloy, whose main constituents are cobalt, chromium, tungsten, and nickel may suffer an allergic response to this implant.

- Once fully deployed, the stent cannot be repositioned.

- Stenting across a bifurcation could compromise future diagnostic or therapeutic procedures.

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The safety and efficacy of the PALMAZ BLUE .018 Peripheral Stent System have not been established for use in patients for treatment of atherosclerotic disease of the femoral artery.

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- When stenting renal arteries, exercise great care to reduce the risk of plaque embolization.

- Avoid stent placement which would obstruct access to a vital side branch.

- It is recommended that stents not be implanted in patients with confirmed pregnancies.

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• To reduce the potential for vessel damage, the inflated diameter of the balloon should approximate the diameter of the vessel just proximal and distal to the stenosis.

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- The device should only be used by physicians who are trained in interventional techniques such as PTA, placement of stents and transhepatic access.
- Prior to stenting, the PALMAZ BLUE.018 Peripheral Stent System should be examined to verify functionality and integrity.
- The delivery system is not designed for use with power injection systems.
- Recrossing a partially or fully deployed stent with adjunct devices must be performed with extreme caution.
- To avoid the possibility of dissimilar metal corrosion, do not implant stents of different metals in tandem where overlap or contact is possible, with an exception: stents made of cobalt-chromium alloy conforming to ASTM F90 (L605) are compatible with stents made of 316L stainless steel and nickel titanium alloy.
- Do not attempt to remove or readjust the stent on the delivery system.
- To assure full expansion, inflate to at least the recommended nominal pressure as shown on the label.
- Do not exceed the rated burst pressure recommended on the label. The rated burst pressure is based on the results of *in-vitro* testing. At least 99.9% of the balloons (with a 95% confidence), will not burst at or below their rated burst pressure. Use of a pressure monitoring device is recommended to prevent over

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Do not advance or retract the catheter unless the balloon is fully deflated under vacuum.

- The minimal acceptable sheath French size is printed on the package label. Do not attempt to pass the delivery system through a smaller size sheath introducer than indicated on the label.

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- Not intended for precise arterial blood

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pressure monitoring.

- Consider the use of systemic heparinization. Flush all devices entering vascular system with sterile heparinized saline or similar isotonic solution.

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Descriptions or specifications in Cordis printed matter, including this publication, are meant solely to generally describe the product at the time of manufacture and do not constitute any express warranties. Cordis Corporation will not be responsible for

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any direct, incidental, or consequential damages resulting from reuse of the product.

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Instruction for Use Stents PALMAZ BLUE

The medical device is the PALMAZ BLUE on .014" Peripheral Stent System

CONTENTS: One (1) PALMAZ BLUE .014 Peripheral Stent System, one (1) flushing needle, one (1) metal introducer tube, one (1) compliance chart, one (1) stent implant card and one (1) instructions for use.

For stent and balloon dimensions, see the carton label. For stent inner diameter vs. inflation pressure, refer to the compliance chart.

Stent

The PALMAZ BLUE Peripheral Stent is a balloon-expandable, laser cut stent made from cobalt-chromium tubing. The stent is mounted on a Cordis Percutaneous Transluminal Angioplasty (PTA) balloon dilatation catheter.

Delivery System

The Cordis AVIATOR PLUS Percutaneous Transluminal Angioplasty (PTA) Balloon Dilatation Catheter, a catheter with a distal DURALYN inflatable balloon, is used as a delivery system for the PALMAZ BLUE Peripheral Stent. The catheter utilizes a Rapid Exchange design, consisting of a single inflation lumen and a distal guidewire lumen. The guidewire lumen begins at the distal tip and terminates at the guidewire exit port. The guidewire exit port (hole) is approximately 25 cm from the distal tip. The maximum guidewire diameter that may be used is printed on the package label. The catheter tip is tapered to facilitate crossing of tight lesions. The proximal hub is

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used as a balloon inflation port. The balloon is inflated by injecting diluted contrast medium through this hub. Two radiopaque marker bands within the balloon indicate the dilating section of the balloon and aid in stent placement. The nominal balloon diameter and length are printed on the identification band near the hub. The working pressure range for the balloon is between the nominal pressure and the rated burst pressure. All balloons distend to sizes above the nominal size at pressures greater than the nominal pressure. Consult the compliance table incorporated with the product for the stent's inner diameter at a given pressure. The design of this catheter does not incorporate a lumen for distal dye injection. The catheter with a usable length of 142 cm contains two proximal shaft markers located 90 cm and 100 cm from the distal tip. Both markers indicate the relative position of the catheter tip to the distal end of the guiding catheter. An additional marker is located at the distal port exit notch and aids in locating the guidewire exit port

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The Cordis PALMAZ BLUE .014 Peripheral Stent System is intended for use in the treatment of atherosclerotic disease of peripheral arteries below the aortic arch.

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The Cordis PALMAZ BLUE .014" Peripheral Stent System is not indicated for use in the arcus aorta or the aorta descendens to the bifurcation aortae.

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- Patients with highly calcified lesions resistant to PTA.
- Patients with a target lesion with a large amount of adjacent acute or sub-acute thrombus.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy.
- Patients with perforated vessels evidenced by extravasation of contrast media.
- A lesion that is within an aneurysm or an aneurysm with a proximal or distal segment adjacent to the lesion.

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- Aneurysm
- Allergic reaction to stent material cobalt chromium or one of its components, or to contrast medium
- Arteriovenous fistula
- Cerebrovascular accident/Stroke
- Death
- Embolization (air, tissue, stent, atherosclerotic or thrombotic material, cholesterol)
- Emergency surgery to correct vascular complications (revascularization/ bypass)
- Failure to deliver the stent to the intended

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а. The patient may be started on anticoagulation medication prior to the procedure if deemed appropriate by the physician.

б. The percutaneous placement of the stent should be done in an angiography procedure room. Angiography should be performed to map out the extent of the lesion. If thrombus is present or suspected, thrombolysis should precede stent deployment using standard practice. Access vessel must be sufficiently patent, or sufficiently recanalized, to proceed with further intervention. Patient preparation and sterile precautions should be the same as for any angioplasty procedure.

2. Selection of Stent Size

а. Measure the length of the target lesion to determine the length of stent required. Size the stent length to extend slightly proximal and distal to the lesion, taking stent foreshortening into account (refer to label for stent length when fully expanded).

б. The appropriate stent length should be selected based on covering the entire obstructed segment with a single stent. **Note:** Should more than one stent be required, place the stent most distal from the puncture site first, followed by placement of the proximal stent in tandem.

с. Measure the diameter of the reference vessel to determine the appropriate size stent system.

3. Preparation of Stent System a. Remove the inner package from the box.

Open the inner package and carefully extract the packaging tube with the stent system and the packaging tray containing flushing needle and introducer tube.

с. Hold the packaging tube in one hand. With the other hand, gently grasp the hub and carefully remove the stent system from the tube.

д. Remove the introducer tube, flushing needle and compliance chart from the tray;

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discard the tray and packaging tube.

е. Inspect the crimped stent for adherence to the balloon. Do not reposition the stent or hand crimp.

ф. Attach a syringe filled with sterile heparinized saline or similar isotonic solution to the flushing needle. Remove the protection sheath from the flushing needle, insert the needle into the tip of the catheter and flush the guidewire lumen.

Caution: Avoid manipulation of stent during flushing of guidewire lumen, as this may disrupt the placement of the stent on the balloon.

г. Attach a three-way stopcock to the catheter's inflation port.

h. Purge the air from a partially filled syringe and connect the syringe to the stopcock.

i. Open the stopcock and induce negative pressure.

j. Hold the syringe and proximal end of the catheter vertically with the balloon tip pointing down.

k. While maintaining the negative pressure, close the stopcock to the inflation port.

l. Remove the syringe.

m. To ensure all air is removed from the balloon and inflation lumen, repeat steps h-l.

n. Prepare an inflation device with diluted contrast medium.

o. Purge the air from the inflation device and connect the inflation device to the stopcock that is connected to the catheter inflation port.

Caution: Non-ionic contrast medium has higher viscosity and precipitation levels than does the ionic type, which may prolong inflation/deflation times.

p. Open the stopcock to the catheter, the inflation lumen and the balloon will slowly be filled with diluted contrast medium.

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Caution: Do not apply negative or positive pressure to the balloon at this time.

Note: After prepping, the catheter can be kept in a coiled configuration by looping the catheter and inserting the proximal shaft into the hub clip.

Caution: Do not clip the distal section of the catheter.

4. Insertion of Guidewire Through Catheter Sheath Introducer (CSI) and Guiding Catheter (GC)

a. Access should be obtained in standard fashion through the femoral approach. The brachial route may be used in cases where the anatomy dictates this. Gain access at the appropriate site utilizing the CSI of the appropriate size allowing for later GC introduction.

Caution: Always use a Cordis CSI or Cordis GC for the implant procedure to protect the puncture site and to avoid dislodging the stent from the balloon.

b. Prepare the guiding catheter recommended on the label.

c. Connect a Tuohy-Borst type hemostasis device to the guiding catheter.

d. Prior to use, flush the guiding catheter lumen with a heparinized saline solution.

e. Introduce the CSI and guiding catheter into the vasculature using the introduction technique of choice.

f. Introduce the guidewire. Advance the guidewire across the target lesion. Close the hemostasis valve around the guidewire. Maintain lesion access with the guidewire across the lesion and with the GC positioned near the target lesion.

Caution: Use care during advancement of the guiding catheter and guidewire to prevent bleeding, dissection or patient discomfort.

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Note: Compatible guidewire diameter is listed on the label.

Note: If the physician determines that predilatation is necessary, standard PTA techniques may be used.

Caution: If PTA technique is used, avoid over expansion of the balloon to prevent bleeding, dissection or patient discomfort.

5. Introduction of Stent System

a. Backload the stent system onto the guidewire. The delivery system has a rapid exchange design. During backloading, the guidewire will exit the stent system at the guidewire exit port. Advance the tip of the stent system to the hemostasis device of the guiding catheter.

b. Open the hemostasis valve as wide as possible and carefully advance the stent system over the guidewire until the balloon section of the stent system is completely introduced into the guiding catheter. Look for and confirm back flow over the balloon. Adjust the hemostasis valve to maintain a snug seal.

c. Advance the stent system until the guidewire exit port has passed the hemostasis valve. Again adjust the hemostasis valve to maintain a snug seal. Continue to advance the stent system over the guidewire to the end of the guiding catheter, while maintaining guidewire position across the lesion.

Caution: When there is a tight fit between the balloon section of the stent system and the guiding catheter, a failure to maintain a snug seal of the hemostasis valve over the stent system during stent system insertion and advancement, could result in air introduction and air entrainment. **Note:** In case a long CSI is used instead of a GC, use the CSI size recommended on the label; insert the stainless steel introducer tube included in the packaging through the hemostasis valve of the CSI. This tube facilitates introduction of the stent/balloon assembly into the CSI and prevents the hemostasis valve from potentially damaging or stripping the

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стент from the balloon. Remove the introducer tube when the stent/ balloon assembly has completely passed the valve (i.e., has advanced to the body of the CSI).

6. Stent Deployment

a. Keeping the GC immobile, observe fluoroscopically as the stent system is advanced through the GC over the guidewire to the target lesion. Carefully cross the lesion with the stent. **Caution:** If strong resistance is met during advancement or withdrawal of the stent system, discontinue movement and determine the cause of resistance before proceeding. If the cause of resistance cannot be determined, secure the stent system to the GC; then remove the GC and stent system as a single unit.

b. Under fluoroscopy, use the balloon marker bands and the radiopaque stent to properly position the stent within the lesion. During positioning, verify that the stent is still centered within the balloon marker bands and has not dislodged.

c. Prior to stent expansion, ensure that stent and balloon are completely exposed from the GC. **Caution:** Never advance the GC over an exposed stent to avoid dislodging the stent. **Caution:** Avoid over tightening of the hemostasis valve since this may restrict the flow of contrast media into and out of the balloon, thereby prolonging inflation/deflation times.

d. Using the inflation device, steadily inflate the balloon under fluoroscopy to the nominal pressure recommended on the label and compliance chart. Expand the diameter of the stent to the diameter of the reference vessel.

Note: Consult the compliance chart for the expanded stent inner diameter at various

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inflation pressures.

Caution: Do not exceed the rated burst pressure shown on the label and compliance chart.

Note: Never use air or any gaseous medium to inflate the balloon.

Note: It is strongly recommended that the guidewire, the stent system, or both, remain across the lesion until the deployment procedure is completed.

7. Delivery System Withdrawal

a. After deploying the stent, while observing under fluoroscopy deflate the balloon by pulling negative pressure with the inflation device, allowing adequate time for the balloon to fully deflate prior to removal.

b. While maintaining negative pressure on the balloon, slightly open the hemostasis valve and carefully withdraw the balloon from the stent. Observe removal of the balloon under fluoroscopy to ensure that the balloon disengages from the stent.

c. Open the hemostasis valve as wide as possible and carefully withdraw and remove the delivery system from the GC while keeping the guidewire in place.

d. Close the hemostasis valve to maintain a snug seal around the guidewire.

e. A post-stent angiogram should be obtained. After a satisfactory post-stent angiogram is obtained, remove the guiding catheter and the guidewire and finally the CSI from the body. Discard the devices.

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Note: Packaging disposal can be handled per normal industrial standard. The used devices shall be handled according to hospital procedures.

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Magnetic Resonance Imaging (MRI): Through non-clinical testing, the Cordis PALMAZ BLUE Peripheral Stent has been shown to be MRI conditional, when tested either individually or in combination of 2 stents with an overlap of 4 mm, at field strengths of 3.0 Tesla or less, a maximum gradient of 5 T/m, and a maximum whole body averaged specific absorption rate (SAR) of 4 W/Kg for 15 minutes of MRI. The PALMAZ BLUE Peripheral Stent should not migrate in this MRI environment. MRI at 3 Tesla or less may be performed immediately following the implantation of the stent. Non-clinical testing has not been performed to rule out the possibility of stent migration at field strengths higher than 3 Tesla. In this non-clinical testing, the stents produced a temperature rise of less than 1.0°C at a maximum whole body averaged specific absorption rate (SAR) of 4 W/Kg for 15 minutes of MRI. The effect of heating in the MRI environment for fractured struts is not known. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the stent. Optimization of MR imaging parameters is recommended.

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The device should only be used by physicians who are trained in interventional techniques and placement of stents.

- Before and during the procedure, appropriate antiplatelet and anticoagulant therapy must be provided to the patient as needed. The therapy should be continued after the procedure for a period of time as determined by the physician.
- Patients with uncorrected bleeding disorders or patients who cannot receive anticoagulation or antiplatelet aggregation therapy may not be suitable for this implant.
- In patients requiring the use of antacids and/or H₂ antagonists before or immediately after stent placement, oral absorption of antiplatelet agents may be adversely affected.
- It is recommended that stents not be implanted in patients with confirmed pregnancies.
- Patients with lesions in arteries to transplanted or bypassed kidneys are not suitable candidates for stent implantation.
- Patients with allergic reactions to a cobalt-chromium alloy, whose main constituents are cobalt, chromium, tungsten, and nickel, may suffer an allergic response to this implant.
- Patients with highly calcified arterial lesions highly resistant to PTA may not be suitable for this implant.
- Caution should be taken with patients with poor renal function who, in the physician's opinion, may be at risk for a contrast medium reaction.
- Patients whose target lesion contains a large amount of acute or subacute thrombus may not be suitable for this implant.

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тween the balloon section of the balloon catheter and the sheath or guiding catheter may cause a risk for introduction of air and air entrapment during advancement of the balloon catheter through the sheath or guiding catheter.

- Stenting across a bifurcation could compromise future diagnostic or therapeutic procedures.
- The safety and efficacy of this product has not been established for use in patients for treatment of atherosclerotic disease of the femoral artery.
- Avoid stent placement that would obstruct access to a vital side branch.

The use of this product for procedures other than those indicated in these instructions is not recommended.

- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.
- Store in a cool, dark and dry place.
- Use the product prior to the "Use By" date specified on the package.

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be manipulated only under fluoroscopy. Do not advance or retract the catheter unless the balloon is fully deflated under vacuum. If resistance is met during manipulation, determine the cause of the resistance before proceeding.

- To assure full expansion, inflate to at least the recommended nominal pressure as shown on the label and compliance chart.
- Do not exceed the rated burst pressure recommended on the label and compliance chart. The rated burst pressure is based on the results of *in-vitro* testing. At least 99.9% of the balloons (with a 95% confidence) will not burst at or below their rated burst pressure. Use of a pressure-monitoring device is recommended to prevent over pressurization and potential complications.
- When treating multiple lesions, the lesion distal to the puncture site should be initially stented, followed by stenting of the proximal lesion. Stenting in this order obviates the need to cross the proximal stent reducing the chance for stent dislodgment.
- As with any type of intravascular implant, infection, secondary to contamination of the stent, may lead to thrombosis, pseudoaneurysm or rupture into a neighboring organ or retroperitoneum.
- In the event of complications, surgical removal of the stent may be required. Standard surgical procedure is appropriate.

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Cordis medical devices may be transported using all available modes of transportation in covered vehicles, according to the transportation rules applicable to each transportation mode, at temperatures ranging from -40°C to $+50^{\circ}\text{C}$ and at 20% - 80% relative humidity.

Cordis medical devices should be stored at temperatures ranging from $+15^{\circ}\text{C}$ to $+25^{\circ}\text{C}$, in a place protected from direct sunlight. Cordis medical devices being in operation are stable in temperature conditions from $+32^{\circ}\text{C}$ to $+42^{\circ}\text{C}$ and to biological liquids.

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The products which had a contact with blood and / or body fluids are subject to disposal and / or destruction according to SanPiN 2.1.7.2790-10 (B class).

Unused products (which had no contact with blood and / or body fluids), including those with expired shelf-life, are subject to disposal and / or destruction, according to SanPiN 2.1.7.2790-10 (Class A).

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Instruction for Use Stents

PALMAZ Peripheral

The PALMAZ balloon expandable Stent for *Large Vessels* is a 316L slotted stainless steel tube that is used to maintain patency in the vessels. The stent has an outside diameter of 4.5 mm uncrimped and is available in two lengths, 40mm and 50mm. When crimped onto a balloon catheter, the stent has a nominal outside diameter of 4.0mm and can be expanded to between 14.0 and 25.0mm in diameter at a balloon inflation pressure of 1.5 to 2 atmospheres. The final diameter of the expanded stent is determined by the size of the fully inflated balloon.

The PALMAZ balloon-expandable Stent for *Large Vessels* is recommended for use with the Mansfield™ Valvuloplasty 15 mm Balloon Catheter (Product Code 2614) for initial deployment of the stent. When the stent is crimped onto the 15mm balloon, the outside diameter of the stent becomes 3.8 mm. If further expansion of the stent is necessary, the Mansfield Valvuloplasty 20 mm Balloon Catheter (Product Code 2616) may be used.

EXPANDED STENT			
P4014		P5014	
Diameter mm	Length mm	Diameter mm	Length mm
15	37.4	15	46.5
20	34.7	20	43.2

After the *Large Vessel* stent has been expanded and properly implanted, the atherosclerotic plaque is compressed and held against the stretched media and adventitia. The stainless steel mesh may be incorporated into the arterial wall as the endothelium bridges over the struts and completely covers the inner stent surface. Experimental evidence has shown that similar metal mesh stents are completely covered with endothelial cells in one week in rabbits and three weeks in dogs.

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The PALMAZ balloon-expandable Stent is indicated for use in peripheral vessels below the aortic arch in patients who are acceptable candidates for percutaneous transluminal angioplasty. The lumen diameter of all vessels between the entry site and the lesion site must be large enough to accept a 12F or 13F introducer sheath. If the vessels are not large enough due to a lesion, an angioplasty procedure must be performed to produce a lumen of sufficient size. The target vessel must be large enough to accommodate the stent expanded to 15mm. In order to perform a successful stent implantation, a CT scan is recommended to ensure the accurate measurement of the passage vessels and the target vessel.

The PALMAZ balloon-expandable stent is contraindicated for use in patients who have the following conditions:

1. Total occlusions of both iliac arteries when stenting via the trans-femoral approach.
2. Marked tortuosity of the iliac artery which prevents passage of the stent when stenting via the trans-femoral approach.
3. Vessel Lumens too small to accommodate introduction of a 12F or 13F introducer sheath.
4. An aneurysmal dilatation of the aorta proximal or distal to the target

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As with any type of vascular implant, infection secondary to contamination of the stent may lead to thrombosis, pseudoaneurysm, rupture into a neighboring organ or the retroperitoneum. The stainless steel stent may cause a thrombus, distal embolization, or may migrate from the site of implant down the arterial lumen. Overstretching of the artery may result in rupture and life threatening bleeding. In addition, since the stent is made of 316L stainless steel, persons with allergic reactions to 316L stainless steel may suffer an allergic response to this implant.

This medical device is used only in the conditions of medical and therapeutic institutions by qualified medical personnel

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1. Initial Angioplasty

- a) After local anesthesia is administered, the arterial or venous system is entered with a puncture needle.
- b) A stiff guidewire is introduced through the needle and is advanced across the stenosis.
- c) The needle is removed and a straight catheter is introduced and advanced over the guidewire.
- d) Inject contrast media through the catheter to confirm the intraluminal position.
- e) While maintaining the guidewire position, exchange the straight catheter for the 12F or 13F introducer sheath.
- f) Select an angioplasty balloon catheter which corresponds to the diameter of the vessel proximal to the lesion. When there is doubt about the distensibility of the lesion, the smallest diameter balloon catheter should be used for the initial dilation.
- g) Connect the side arm adaptor of the introducer sheath to a pressure transducer to record the arterial pressure distal to the obstruction. An initial dilation of the lesion should be made with an appropriately sized balloon catheter.
- h) Following dilatation of the lesion, an arteriographic image should be recorded (1) to determine the adequacy of the dilation and (2) to place an external metal marker at the lesion site.

2. Stent Preparation

- a) Use sterile technique during preparation and deployment of the stent. The PALMAZ balloon-expandable Stent is provided STERILE for one use only. The stainless steel Introducer Tube is provided NON-STERILE and must be sterilized prior to use. The Introducer Tube may be sterilized and resterilized by autoclaving in accordance with hospital procedures.
- b) Use the recommended 15mm balloon catheter for initial stent deployment.

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Only use the 20mm balloon catheter if post-implantation expansion is necessary.

- c) After careful inspection of the stent package looking for damage to the sterile barrier, remove the stent from the package and rinse in sterile heparinized saline.
- d) Remove the 15mm balloon catheter (Mansfield Product Code 2614) from the sterile package in preparation for stent mounting.
- e) Following the catheter manufacturer's guidelines, prepare the balloon catheter. Inflate the balloon with a mixture of saline solution and contrast media.
- f) Deflate the balloon by pulling negative pressure. While holding negative pressure, close the stopcock. Maintain negative pressure until the stent is mounted.
- g) Carefully refold the balloon to a low profile. Grasp balloon wings and pull gently to center the catheter shaft and allow for the smallest balloon profile.
- h) Slide the stent over the distal end of the balloon and center on the balloon. The fully inflated balloon length is 55mm; the stent lengths are 40mm and 50mm. If an encounter with a bifurcation near the implant site is expected, the stent may be slightly adjusted proximally on the balloon. Care must be taken not to overadjust the stent. An overadjustment could cause the misshapen expansion of the stent or cause its slippage off the balloon.
- i) Insert a stiff guidewire or mandrel in the lumen of the balloon for protection during crimping.
- j) Manually finger crimp the stent onto the balloon. Do not roll the stent on the balloon or turn the balloon while crimping. Firmly crimp the stent in one location and then rotate your fingers a 1/8 turn around it and crimp again. Continue to firmly crimp the stent rotating a 1/8 turn after each crimp until the stent struts have a uniform "bow tie" shape. This will assure that the stent is securely crimped to the balloon.
- k) Release the negative pressure and remove the stiff guidewire.
- l) Gently tug the stent to be sure it is securely crimped. If the stent feels loose, repeat step j. Repeat test.

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m) Pull negative pressure to check the balloon. If air bubbles appear, repeat steps d–n with a new balloon catheter. If no bubbles appear, slowly release the negative pressure.

n) Gently tug the stent again to be sure it is securely crimped. If the stent feels loose, use your fingers to further crimp it.

3. Introduction of the Stent

a) With the introducer sheath still in place, carefully slide the flared end of the metal introducer tube (INTR1213F) over the mounted stent. Insert the beveled end of the metal tube into the hemostatic valve of the introducer sheath until it stops. Slide the stent-mounted catheter through the tube until the stent is completely into the introducer sheath. Slide the metal tube back to the proximal end of the catheter and close the hemostatic valve. The metal tube can remain at the proximal end of the catheter shaft.

b) Under fluoroscopy, advance the stent-mounted balloon assembly to the dilated lesion using the external metal marker as a guide.

c) With the stent in position across the lesion, withdraw the introducer sheath to fully expose the stent and the complete length of the balloon.

d) Inject a small amount of contrast media into the introducer sheath to assess the correct positioning. (Note: Complications could arise if the introducer sheath is not long enough to cross the lesion or if the site of entry is distal to the lesion. In this

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case, other means of assessing correct positioning of the stent must be employed.)

е) Deploy the stent by slowly inflating the balloon with a mixture of saline solution and contrast media until the balloon appears dumbbell shaped. The inflation speed may now be increased. The stent is designed to expand fully at inflation pressures of 1.5 to 2 atmospheres. Complete expansion of the stent should be determined by fluoroscopy. Do not exceed the catheter manufacturer's recommendations for maximum balloon inflation pressure.

ф) Deflate the balloon.

г) Rotate and advance slightly the deflated balloon delivery catheter to ensure the stent is free and properly deployed.

h) Remove the delivery catheter.

NOTE: If the patient has significant femoral disease or if additional surgical intervention is contemplated, the stent may be placed in an operating room following surgical cut-down and isolation of the femoral artery.

NOTE: The diameter of the stent may be increased post-placement by expanding with a large diameter balloon. Do not exceed the maximum recommended expanded stent diameter of 25.0 mm.

NOTE: The stent is adequate for use in lesions of approximately 3.5 or 4.5 cm in length, exact length depends on final diameter. In longer lesions, a second stent may be placed slightly overlapping the previous one. Due to the small number of multiple stent implants reports to date, no judgement can be made concerning the relative efficacy of single and multiple stent implants. Excessive angulation between overlapped stents due to tortuosity in the aorta should be avoided.

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4. Post Stent Placement Procedure

a) A post-stenting angiogram should be taken by placing a 6.5F straight catheter through the introducer sheath using a guidewire as a leader to avoid penetration of the vessel and stent dislodgment.

b) Remove the catheter, guidewire, and introducer sheath. Attain hemostasis using standard practices.

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The PALMAZ balloon-expandable Stent is supplied sterile for one use only. Use the PALMAZ balloon-expandable Stent prior to the "Use By" date printed on the package. Use only the recommended Boston Scientific Mansfield Valvuloplasty catheter as a stent delivery catheter. Maximum balloon inflation pressure must not exceed the recommended inflation pressure specified in the manufacturer's instructions. Use of an inflation device with a manometer is recommended during this procedure so as to minimize the risk of overinflation. The inflation diameter of the balloon used during stent delivery should approximate the diameter of the dilated lesion and intended arterial implant site. To assure full expansion, the balloon should be inflated to at least 1.5 to 2 atmospheres of pressure. The stent and delivery balloon should be used only by physicians trained in PTA therapy and in stent implantation.

Antiplatelet or anticoagulant therapy may be used as deemed appropriate by the physician. When catheters are in the body, they should be manipulated only under fluoroscopy. Radiographic equipment that provides high quality images is needed. Circumferential tear of the carrier balloon distal to the stent and prior to complete expansion of the stent may cause the balloon to

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become tethered to the stent, requiring surgical removal. In case of rupture of an adequately sized balloon, it can be withdrawn and a new balloon catheter exchanged over a guidewire to complete expansion of the stent. The fate of branch vessel which are crossed by the stent is unknown at present. The long-term outcome following repeat dilatation of endothelialized stents is unknown at present. The stent may cause artifacts with MRI scans due to distortion of the magnetic field. The artifacts caused by the 316L stainless steel stent should not be greater than that caused by surgical metal clips. A MRI scan should not be used until the implanted stent has had a chance to heal (estimated to be 8 weeks), in order to minimize the risk of migration of the stent under a strong magnetic field.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s)

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described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

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Cordis medical devices may be transported using all available modes of transportation in covered vehicles, according to the transportation rules applicable to each transportation mode, at temperatures ranging from - 40 °C to +50 °C and at 20% - 80% relative humidity.

Cordis medical devices should be stored at temperatures ranging from +15°C to + 25°C, in a place protected from direct sunlight. Cordis medical devices being in operation are stable in temperature conditions from +32 °C to +42 °C and to biological liquids.

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The Cordis TRAPEASE Permanent Vena Cava Filter and Introduction Kit is designed for percutaneous delivery of a permanent vena cava filter to the inferior vena cava (IVC).

The symmetrical and self-centering TRAPEASE Filter is laser cut from nickel titanium alloy (Nitinol) tubing. The proximal and distal baskets of TRAPEASE Filter, which consist of struts in a six diamond shape configuration, are designed for optimal clot capture. The baskets are connected by six straight struts which contain a proximal and distal hook designed for fixation of TRAPEASE Filter to the vessel wall (Figure 1).

The constrained filter is flexible and achieves its unconstrained diameter upon deployment in the IVC. Upon deployment, the filter imparts an outward radial force on the luminal surface of the vena cava to ensure proper positioning and stability. TRAPEASE Filter is designed to prevent pulmonary embolism while maintaining caval patency.

The introduction kit is comprised of a Cordis 6F (2.0 mm) BRITE TIP Catheter Sheath Introducer (Sheath Introducer), an obturator with marker and a 0.035" (0.89 mm) J-guidewire. The obturator serves to advance the filter through the Sheath Introducer to the implantation site. The constrained filter is supplied in a plastic storage tube which is to be loaded as a system into the Sheath Introducer hemostasis valve.

The length of the introduction kit and the recommended venous access sites are given in Table I.

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TRAPEASE Filter is indicated for the prevention of recurrent pulmonary embolism via percutaneous placement in the inferior vena cava in the following situations:

- Pulmonary thromboembolism when anticoagulants are contraindicated.
- Failure of anticoagulant therapy in thromboembolic diseases.
- Emergency treatment following massive pulmonary embolism where anticipated benefits of conventional therapy are reduced.
- Chronic, recurrent pulmonary embolism where anticoagulant therapy has failed or is contraindicated.

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- Patients with risk of septic embolism.
- Patients with uncontrolled infectious disease.
- Patients with an IVC diameter larger than 30 mm.
- Patients contraindicated for procedures under fluoroscopy.
- Patients with demonstrated hypersensitivity to one of the components of TRAPEASE Filter

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Procedures requiring percutaneous interventional techniques should not be attempted by physicians unfamiliar with the possible complications.

Complications may occur at any time during the procedure.

Possible procedure complications include, but are not limited to, the following:

- air embolism
- hematoma at the puncture site
- incorrect positioning of the filter
- perforation of the vessel wall
- restriction of blood flow
- occlusion of small vessel
- distal embolization
- infection
- intimal tear
- filter fracture
- thrombus formation.

Possible longer-term complications associated with filter implantation include, but are not limited to, the following:

- filter obstruction
- filter perforation of the vena cava wall
- filter migration
- filter fracture
- recurrent pulmonary embolism.

Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut.

However, reports of adverse clinical sequelae from filter fractures are rare.

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formation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.

- TRAPEASE Filter should only be used by physicians who are trained in percutaneous interventional techniques, for instance placement of vena cava filters. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.

- Persons with allergic reactions to nickel titanium (nitinol) may suffer an allergic response to this implant.

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- Store in a cool, dark, dry place.

- Do not use if package is open or damaged.

- Use prior to "Use By" date.

- Do not autoclave or resterilize by any other method. Do not expose to temperatures above 54°C (130°F) for any length of time.

- Do not expose to organic solvents.

- TRAPEASE Filter has been tested and qualified with the accompanying accessories. The use of any other accessory could result in complications and/or an unsuccessful procedure. If strong resistance is met during any stage of the procedure, discontinue the procedure and determine the cause before proceeding.






- The TRAPEASE Filter has been shown to be MRI safe at field strengths of 3.0 Tesla or less, a maximum gradient of 5 T/m, and a maximum whole body averaged specific absorption rate (SAR) of 4 W/Kg for 15 minutes of MRI through non-clinical testing. The TRAPEASE Filter should not migrate in this MRI environment. MRI at 3

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Tesla or less may be performed immediately following the implantation of the filter. Non-clinical testing has not been performed to rule out the possibility of filter migration at field strengths higher than 3 Tesla. In testing, filters produced a temperature rise of less than 1.8°C at a maximum whole body averaged specific absorption rate (SAR) of 4 W/Kg for 15 minutes of MRI. MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the filter.

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1. Select a suitable venous access site.

Note: It is preferable to use the right-sided approach because filter advance can be problematic when using a tortuous left-sided vessel approach.

2. Remove the components of the introduction kit from the package using sterile technique.

3. Flush the guidewire with sterile heparinized saline or suitable isotonic solution via a syringe connected to the luer hub of the guidewire dispenser.

4. Remove air from the Sheath Introducer by flushing with heparinized saline or suitable isotonic solution.

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5. Insert the vessel dilator through the Sheath Introducer hemostasis valve, snapping it into place at the hub (Figure 3). Flush with heparinized saline or suitable isotonic solution.

6. Introduce the cannula of an angiographic needle into the vessel using aseptic technique. To straighten the J tip of the guidewire for insertion, slide the guidewire introducer over the J tip. Holding the needle in place, insert the guidewire through the needle and into the vessel. Gently advance the guidewire to the desired location.

Caution: Do not withdraw a PTFE-coated guidewire through a metal cannula. Withdrawing may damage the guidewire coating.

7. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Sheath Introducer is inserted into the vasculature. Pass the Sheath Introducer over the guidewire, grasping the Sheath Introducer close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue and into the vessel (Figure 4).

8. Position the tip of the Sheath Introducer in the IVC below renal veins (the intended implantation site). Locate the BRITE TIP (radiopaque) tip and/or use contrast medium to determine the position of the sheath introducer tip in the vasculature.

Caution: Do not use with Ethiodol or Lipiodol contrast media, or other such contrast media, which incorporate the components of these agents.

9. Detach the vessel dilator from the Sheath Introducer by releasing the snap-fit ring at the hub (Figure 5). Withdraw the guidewire and dilator.

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Caution: Do not leave a Sheath Introducer in place for extended periods of time without a dilator or an obturator to support the

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cannula wall.

10. Aspirate from the sideport extension to remove any potential air.

11. Determine the diameter of the vena cava at the intended implantation site. The filter is intended for an IVC of ≤ 30 mm.

12. Place the storage tube (containing filter) as far as possible into the Sheath Introducer hemostasis valve (Figure 6)

13. Slowly advance the filter into the Sheath Introducer by advancing the obturator until the filter is positioned well into the cannula of the Sheath Introducer.

14. To restore hemostasis, withdraw the storage tube out of the Sheath Introducer hemostasis valve. The storage tube can remain on the obturator.

15. Continue to advance the filter until the marker on the obturator is positioned at the Sheath Introducer hemostasis valve. This indicates that the filter is at the distal tip of the Sheath Introducer but still fully in the Sheath Introducer (Figure 7).

Note: If filter advance is problematic when using a tortuous vessel approach, stop filter advance prior to the curve. Advance the Sheath Introducer to negotiate the curve, then continue to advance the filter.

16. Filter release is performed under continuous fluoroscopy. Prior to releasing the filter from the Sheath Introducer ensure that the intended filter location in the IVC is correct. Should the filter location be incorrect, reposition the Sheath Introducer.

Note: The constrained, unexpanded, length of the filter is 64 mm. During release from the Sheath Introducer, the length of the filter will shorten as the filter expands. The ends of the filter will shorten equally as they move toward the axial center of the filter. Filter foreshortening is depicted in Figure 2. The length of the deployed filter is deter-

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mined by its expanded diameter in the IVC. The shortening at each end of the expanded filter and the corresponding filter length is provided in Table II.

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17. Keeping the obturator fixed in position, pull back the Sheath Introducer over the obturator to release the filter from the Sheath Introducer (Figure 8).

18. Ensure that the filter is fully released.

19. Once the filter has been released, remove the obturator from the Sheath Introducer.

20. Perform a control cavography before terminating the procedure.

21. Remove the Sheath Introducer when clinically indicated by placing compression on the vessel above the puncture site, and slowly withdraw the Sheath Introducer.

22. Discard the introduction kit and packaging materials. **Note:** After use, the introduction kit and packaging materials may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable local, state and federal laws and regulations.

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Note: Filter values given are approximate.

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There is no express or implied warranty, including without limitation any implied warranty of merchantability or fitness for a particular purpose, on the Cordis product(s) described in this publication. Under no circumstances shall Cordis be liable for any direct, incidental, or consequential damages

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other than as expressly provided by specific law. No person has the authority to bind Cordis to any representation or warranty except as specifically set forth herein.

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The Cordis OPTEASE Retrieable Vena Cava Filter (Retrieable Filter) is designed for percutaneous delivery of a retrievable vena cava filter to the inferior vena cava (IVC).

The OPTEASE Retrieable Filter can be retrieved within a specified period after implantation or remain implanted as a permanent filter.

The self-centering OPTEASE Retrieable Filter is laser cut from nickel titanium alloy (Nitinol) tubing. The proximal and distal baskets of the OPTEASE Retrieable Filter, which consist of struts in a six diamond shape configuration, are designed for optimal clot capture. The baskets are connected by six straight struts. A single row of fixation barbs is present at the cranial end of the struts. These barbs, intended for fixation to the vessel wall, are extensions of the parallel struts. A hook is centrally located at the caudal basket extremity and allows for filter retrieval using a snare (Figure 1).

The constrained filter is flexible and achieves its unconstrained diameter upon deployment in the IVC. Upon deployment, the filter imparts an outward radial force on the luminal surface of the vena cava to ensure proper positioning and stability. The OPTEASE Retrieable Filter is designed to prevent pulmonary embolism while maintaining caval patency.

The introduction kit is comprised of a filter storage tube, a 6F (2.0 mm ID) Cordis BRITE TIP Catheter Sheath Introducer (Sheath Introducer), an Angiographic Vessel Dilator with an open end, side holes and two marker bands at a distance of 30 mm (end-to-end), and an obturator with marker and storage tube securement bump.

The Angiographic Vessel Dilator is designed to also provide angiographic visualization and linear measurement of the IVC when combined with the delivery of radiopaque contrast media. The obturator serves to advance the filter from the storage tube into the Sheath Introducer and to advance the filter through the Sheath Introducer to the implantation site. The constrained filter is supplied in a plastic storage tube, which is to be loaded as a system into the

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Sheath Introducer hemostasis valve. The storage tube is printed with arrows and text to indicate the access site orientation. The arrows will point into the Sheath Introducer hemostasis valve and can only be inserted in this orientation. The filter must be deployed in the patient with the hook oriented in the caudal position. In this orientation the fixation barbs are designed to prevent the filter from migrating towards the heart and allow retrieval of the filter via the femoral vein.

The length of the introduction kit, the recommended guidewire length and diameter, and the recommended venous access sites are given in Table I.

Retrieval of the OPTEASE Retrievable Filter is from femoral venous access only

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Catalog Number	Recommended Guidewire Length (cm)	Diameter (In.)	Recommended Venous Access Site
466-F210AJ	150	0.035" (0.89 mm) Cat. No. 502-521	Jugular

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The OPTEASER Retrievable Filter is indicated for the prevention of pulmonary embolism (PE) via percutaneous placement in the IVC in patients considered at high risk of PE.

The Angiographic Vessel Dilator is designed to provide angiographic visualization and linear measurement of the vasculature when combined with the delivery of radiopaque contrast media to the vena cava.

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Vena cava filters should not be implanted in:

- Patients with risk of septic embolism.
- Patients with uncontrolled infectious disease.
- Patients with an IVC diameter larger than 30 mm.
- Patients contraindicated for procedures under fluoroscopy.
- Patients with demonstrated hypersensitivity to one of the components of the OPTEASE Retrievable Filter.

There are no known contraindications for use of the Angiographic Vessel Dilator.

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Procedures requiring percutaneous interventional techniques should not be attempted by physicians unfamiliar with the possible complications. Complications may occur at any time during the procedure. Possible procedure complications include, but are not limited to, the following:

- air embolism
- hematoma at the puncture site
- incorrect positioning of the filter
- incorrect orientation of the filter
- perforation of the vessel wall
- restriction of blood flow
- occlusion of small vessel
- distal embolization
- infection
- intimal tear
- filter fracture
- thrombus formation
- incomplete expansion of the filter.

Possible long-term complications associated with filter implantation include, but are not

limited to, the following:

- filter obstruction/thrombosis
- filter perforation of the vena cava wall
- filter migration
- filter fracture
- recurrent pulmonary embolism
- pain
- peripheral edema.

Implant of the OPTEASE Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures, ineffective pulmonary embolism prevention or death.

Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal os-

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tia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare. Retrieval of the Cordis OPTEASE Retrievable Filter with a filter fracture present may result in complications.

This medical device is used only in the conditions of medical and therapeutic institutions by qualified medical personnel.

- Retrieval of the OPTEASE Retrievable Filter is possible only from femoral vein approach. Do not implant the filter with the intention to retrieve using jugular access.
- The OPTEASER Filter must be implanted with the fixation barbs in the cranial direction and retrieval hook in the caudal direction of the patient.
- Implant of the OPTEASE Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures, ineffective pulmonary embolism prevention or death.
- All components in the OPTEASE Introduction Kit with Angiographic Vessel Dilator are for single use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired through cleaning, resterilization, or reuse and may cause adverse patient reactions. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the OPTEASE Introduction Kit.

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- The OPTEASE Retrieval Filter can be retrieved up to and including 12 days after placement. The OPTEASE Retrieval Filter is considered a *permanent* implant if it is not retrieved within the specified time period.

- Retrieval is performed using the Cordis OPTEASE Retrieval Catheter or equivalent, and the recommended accessories (Table III). These devices used for filter retrieval are not included in the OPTEASE Retrieval Filter Introduction Set. The IFU for OPTEASE Retrieval Filter retrieval is contained in the Cordis OPTEASE Retrieval Catheter packaging.

- When injecting contrast medium through the Angiographic Vessel Dilator, do not exceed maximum pressure rating of 800 psi. Ensure that a high-pressure connection line is used.

- After filter implantation, any catheterization procedure requiring passage of a device through the filter may be impeded.

- The OPTEASE Retrieval Filter is supplied constrained in a plastic storage tube indicating the appropriate orientation for jugular approach. Never reload a fully ejected filter into the storage tube as this could result in incorrect orientation for the selected access site. Never reload a (partially) ejected filter into the storage tube as this could affect its shape, and function. Accordingly, Cordis will not be responsible for any direct, incidental or consequential damages resulting from replacement of the OPTEASE Retrieval Filter in the plastic storage tube.

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- Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.

- The OPTEASE Retrieable Filter should only be used by physicians who are trained in diagnostic and percutaneous interventional techniques, for instance placement of vena cava filters. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.

- Persons with allergic reactions to nickel titanium (Nitinol) may suffer an allergic response to this implant.

- Retrieval of the OPTEASE Retrieable Filter with a filter fracture present may result in complications.

- In patients with acute pulmonary embolism in which the OPTEASE Filter is inserted as an alternative to anticoagulation therapy, a conventional course of anticoagulation therapy should be started if the risk of bleeding resolves.

- Where practicable, a retrievable filter should be removed if a conventional course of anticoagulation therapy can be started, or if the underlying reason for indication is resolved.

- Store in a cool, dark, dry place.
- Do not use if package is open or damaged.
- Use prior to "Use By" date.
- Do not autoclave or resterilize by any other method. Do not expose to temperatures above 54°C (130°F) for any length of time.

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- This product is designed and intended for single use. It is not designed to undergo re-processing and re-sterilization after initial use. Reuse of this product, including after re-processing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

- Do not expose to organic solvents.

- The OPTEASE Retrievable Filter has been tested and qualified with the accompanying accessories. The use of any other accessory could result in complications and/or an unsuccessful procedure.

- Advancement, manipulation and withdrawal of the OPTEASE Filter or its accessories should always be performed under fluoroscopic guidance.

- If strong resistance is met during any stage of the procedure, discontinue the procedure and determine the cause before proceeding.

- The safety and effectiveness of the OPTEASE Filter has not been established in patients who are known to be pregnant or who are lactating.

- **MRI Information** Non-clinical testing has demonstrated that the OPTEASE Retrievable Filter is MR Conditional. It can be scanned safely immediately following the implantation of the filter under the following conditions:

- Static magnetic field of 1.5 or 3.0 Tesla.

- Spatial gradient field of 2500 Gauss/cm (25 Tesla/m) or less.

- Maximum whole-body-averaged specific absorption rate (WB-SAR) of 2 W/kg or 4 W/kg for 15 minutes of MR scanning (i.e. per pulse sequence).

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— Normal operating mode of the MR system. The OPTEASE Retrieable Filter should not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of filter migration at field strengths higher than 3 Tesla.

• MRI-Related Heating

— In non-clinical testing, the OPTEASR Retrieable Filter produced a temperature rise of less than 1.8°C at a maximum WB-SAR of 4 W/Kg for 15 minutes of MR scanning.

— In non-clinical testing using computer modeling, the OPTEASE Retrieable Filter produced an estimated in-vivo temperature rise of less than 0.5°C at a WB-SAR of 2 W/kg for 15 minutes of MR scanning in a 3 Tesla / 128 MHz MR scanner and in a 1.5 Tesla / 64 MHz stand-alone body coil.

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3. Remove air from the Sheath Introducer by flushing with heparinized saline or suitable

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isotonic solution.

4. Insert the Angiographic Vessel Dilator through the Sheath Introducer hemostasis valve, snapping it into place at the hub (Figure 3). Flush with heparinized saline or suitable isotonic solution.

5. A Cordis J-tip 0.035" (0.89 mm) EMERALDR Guidewire is recommended to introduce the Sheath Introducer and the Angiographic Vessel Dilator. The length and the maximum diameter of the recommended guidewire are listed in Table I.

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12. Reintroduce the guidewire. Advance the tip of the Sheath Introducer with the

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Angiographic Vessel Dilator in place to the desired location in the IVC.

13. Remove the guidewire. Detach the Angiographic Vessel Dilator from the Sheath Introducer by releasing the snap-fit ring at the hub (Figure 5). Withdraw the Angiographic Vessel Dilator.

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19. Filter release (or deployment) is performed under continuous fluoroscopy. Prior to releasing the OPTEASE Retrievable

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Filter from the Sheath Introducer ensure that:

a. the intended filter location in the inferior vena cava is correct. Should the filter location be incorrect, reposition the Sheath Introducer; and b. the filter retrieval hook is oriented towards the caudal end of the vena cava. Should the filter orientation be incorrect, remove the Sheath Introducer (with filter inside) from the patient and discard the system. Recommence the procedure using a new sterile Sheath Introducer and storage tube with filter.

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clinically indicated by placing compression on the

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Note: After use, the introduction kit and packaging materials may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and with applicable local, state and federal laws and regulations

The OPTEASE Retrieable Filter can be retrieved up to and including 12 days after implantation from femoral access only.

The recommended devices required for the retrieval of the OPTEASE Retrieable Filter are

listed in Table III below.

The IFU for the retrieval of the OPTEASE Retrieable Filter is contained in the packaging of the Cordis OPTEASE Retrieval Catheter.

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Table II. Filter Length per IVC Diameter
Unexpanded Filter Length = 67 mm

IVC Diameter (mm)	Expanded Filter Length (mm)	Filter Shortening at each end (mm)
10	66	0.5
12	65	1
14	64	1.5
16	64	1.5
18	63	2
20	62	2.5
22	61	3
24	59	4
26	57	5
28	56	5.5
30	54	6.5

Note: Filter length and shortening values given are approximate.

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The OPTEASE Retrieval Filter can be retrieved up to and including 12 days after implantation from femoral access only.

The recommended devices required for the retrieval of the OPTEASE Retrieval Filter are listed in Table III below.

The IFU for the retrieval of the OPTEASE Retrieval Filter is contained in the packaging of the Cordis OPTEASE Retrieval Catheter.

Table III. Recommended Devices for the Retrieval of the OPTEASE® Retrieval Filter

Recommended Devices	Dimension	Cat. No.
OPTEASE® Retrieval Catheter with radiopaque tip or equivalent	10F, 80 cm length	466-C210F
Cordis J-tip EMERALD® 0.035" (0.89 mm) Guidewire or equivalent	260 cm length	502-455
Cordis AVANTI®+ Catheter Sheath Introducer or equivalent	10F (3.3 mm), 11 cm length	504-610X
Goose Neck Snare Kit or equivalent	10 mm loop diameter, 4F catheter	GN1000
	15 mm loop diameter, 6F catheter	GN1500
	20 mm loop diameter, 6F catheter	GN2000
	25 mm loop diameter, 6F catheter	GN2500
	30 mm loop diameter, 6F catheter	GN3000
	All 120 cm snare length and 102 cm catheter length	

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The Angiographic Vessel Dilator is designed to also provide angiographic visualization and linear measurement of the IVC when combined with the delivery of radiopaque contrast media. The obturator serves to advance the filter from the storage tube into the Sheath Introducer and to advance the filter through the Sheath Introducer to the implantation site. The constrained filter is supplied in a plastic storage tube, which is to be loaded as a system into the Sheath Introducer hemostasis valve. The stor-

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age tube is printed with arrows and text to indicate the access site orientation. The arrows will point into the Sheath Introducer hemostasis valve and can only be inserted in this orientation. The filter must be deployed in the patient with the hook oriented in the caudal position. In this orientation the fixation barbs are designed to prevent the filter from migrating towards the heart and allow retrieval of the filter via the femoral vein.

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Catalog Number	Recommended Guidewire Length (cm)	Recommended Guidewire Diameter (in.)	Recommended Venous Access Site
466-F210BJ	260	0.035" (0.89 mm) Cat. No. 502-455	Jugular/Antecubital

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- air embolism
- hematoma at the puncture site
- incorrect positioning of the filter
- incorrect orientation of the filter

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- perforation of the vessel wall
- restriction of blood flow
- occlusion of small vessel
- distal embolization
- infection
- intimal tear
- filter fracture
- thrombus formation
- incomplete expansion of the filter.

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Implant of the OPTease Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures, ineffective pulmonary embolism prevention or death.

Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare. Retrieval of the Cordis OPTease Retrieval Filter with a filter fracture present may result in complications.

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- Retrieval of the OPTEASE Retrieable Filter is possible only from femoral vein approach. Do not implant the filter with the intention to retrieve using jugular access.

- The OPTEASER Filter must be implanted with the fixation barbs in the cranial direction and retrieval hook in the caudal direction of the patient.

- Implant of the OPTEASE Filter with the hook oriented in the cranial direction can result in life threatening or serious injury including, but not limited to dissection, vessel perforation, migration of the filter with secondary damage to cardiac structures, ineffective pulmonary embolism prevention or death.

- All components in the OPTEASE Introduction Kit with Angiographic Vessel Dilator are for single use only. Do not resterilize or reuse. Structural integrity and/or function may be impaired through cleaning, resterilization, or reuse and may cause adverse patient reactions. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from reuse of any of the components in the OPTEASE Introduction Kit.

- The OPTEASE Retrieable Filter can be retrieved up to and including 12 days after placement. The OPTEASE Retrieable Fil-

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тер is considered a *permanent* implant if it is not retrieved within the specified time period.

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- Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from

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filter fractures are rare.

- The OPTEASE Retrieable Filter should only be used by physicians who are trained in diagnostic and percutaneous interventional techniques, for instance placement of vena cava filters. Accordingly, Cordis will not be responsible for any direct or consequential damages or expenses resulting from use by untrained personnel.
- Persons with allergic reactions to nickel titanium (Nitinol) may suffer an allergic response to this implant.
- Retrieval of the OPTEASE Retrieable Filter with a filter fracture present may result in complications.
- In patients with acute pulmonary embolism in which the OPTEASE Filter is inserted as an alternative to anticoagulation therapy, a conventional course of anticoagulation therapy should be started if the risk of bleeding resolves.
- Where practicable, a retrievable filter should be removed if a conventional course of anticoagulation therapy can be started, or if the underlying reason for indication is resolved.

- Store in a cool, dark, dry place.
- Do not use if package is open or damaged.
- Use prior to "Use By" date.
- Do not autoclave or resterilize by any other method. Do not expose to temperatures above 54°C (130°F) for any length of time.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

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- Do not expose to organic solvents.
- The OPTEASE Retrievable Filter has been tested and qualified with the accompanying accessories. The use of any other accessory could result in complications and/or an unsuccessful procedure.
- Advancement, manipulation and withdrawal of the OPTEASE Filter or its accessories should always be performed under fluoroscopic guidance.
- If strong resistance is met during any stage of the procedure, discontinue the procedure and determine the cause before proceeding.
- The safety and effectiveness of the OPTEASER Filter has not been established in patients who are known to be pregnant or who are lactating.

• **MRI Information** Non-clinical testing has demonstrated that the OPTEASE Retrievable Filter is MR Conditional. It can be scanned safely immediately following the implantation of the filter under the following conditions:

- Static magnetic field of 1.5 or 3.0 Tesla.
- Spatial gradient field of 2500 Gauss/cm (25 Tesla/m) or less.
- Maximum whole-body-averaged specific absorption rate (WB-SAR) of 2 W/kg or 4 W/kg for 15 minutes of MR scanning (i.e. per pulse sequence).
- Normal operating mode of the MR system. The OPTEASE Retrievable Filter should not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of filter migration at field strengths higher than 3 Tesla.

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• MRI-Related Heating

— In non-clinical testing, the OPTEASR Retrievable Filter produced a temperature rise of less than 1.8°C at a maximum WB-SAR of 4 W/Kg for 15 minutes of MR scanning.

— In non-clinical testing using computer modeling, the OPTEASE Retrievable Filter produced an estimated in-vivo temperature rise of less than 0.5°C at a WB-SAR of 2 W/kg for 15 minutes of MR scanning in a 3 Tesla / 128 MHz MR scanner and in a 1.5 Tesla / 64 MHz stand-alone body coil.

• MRI-Related Artifacts

MR image quality may be compromised if the area of interest is in the exact same area or relatively close to the position of the filter. The image artifact extends approximately 3mm for the spin echo sequence and 7mm for the gradient echo sequence, both inside and outside the device lumen, when scanned in nonclinical testing in a 3 Tesla / 128 MHz MR scanner

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Note: It may be preferable to use the right-sided approach because filter advancement can be problematic when using a tortuous left-sided vessel approach.

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3. Remove air from the Sheath Introducer by flushing with heparinized saline or suitable isotonic solution.

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5. A Cordis J-tip 0.035" (0.89 mm) EMERALDR Guidewire is recommended to introduce the Sheath Introducer and the Angiographic Vessel Dilator. The length and the

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maximum diameter of the recommended guidewire are listed in Table I.

6. Introduce the cannula of an angiographic needle into the vessel using aseptic technique.

7. Slide the guidewire introducer over the J tip of the guidewire to straighten it for insertion. Holding the needle in place, insert the guidewire through the needle and into the vessel. Gently advance the guidewire to the desired location.

8. Holding the guidewire in place, withdraw the needle and apply pressure to the puncture site until the Sheath Introducer with the Angiographic Vessel Dilator is inserted into the vasculature. Pass the Sheath Introducer system over the guidewire, grasping the Sheath Introducer close to the skin to prevent buckling. Using a rotating motion, advance the assembly through the tissue and into the vessel (Figure 4).

9. Position the Sheath Introducer radiopaque tip and the marker bands of the Angiographic Vessel Dilator in the inferior vena cava below the renal veins in preparation for an angiographic overview of the IVC.

10. Remove the guidewire.

11. Determine the diameter of the inferior vena cava at the intended implantation site below the lowest renal vein by injecting contrast medium through the Angiographic Vessel Dilator and using its marker bands as a reference. The distance between the two marker bands, end-to-end, is 30 mm.

Caution: Do not use with Ethiodol or Lipiodol* contrast media, or other such contrast media, which incorporate the components of these agents.

Caution: Do not exceed 800 psi when injecting through a power injector. Ensure that a high-pressure connection line is used.

12. Reintroduce the guidewire. Advance the tip of the Sheath Introducer with the Angiographic Vessel Dilator in place to the desired location in the IVC.

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Caution: To avoid damage to the Sheath Introducer tip, do not withdraw the dilator until the Sheath Introducer tip is at the de-

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a. NOTE: The OPTEASE Filter is indicated for retrieval from the femoral approach only. The retrieval hook should be positioned caudally to allow the anti-migration fixation hooks to engage the caval wall properly.

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a. the intended filter location in the inferior vena cava is correct. Should the filter location be incorrect, reposition the Sheath Introducer; and **b.** the filter retrieval hook is oriented towards the caudal end of the vena cava. Should the filter orientation be incorrect, remove the Sheath Introducer (with filter inside) from the patient and discard

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the system. Recommence the procedure using a new sterile Sheath Introducer and storage tube with filter.

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25. Discard the introduction kit and packaging materials.

Note: After use, the introduction kit and packaging materials may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and

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with applicable local, state and federal laws and regulations

The OPTEASE Retrievable Filter can be retrieved up to and including 12 days after implantation from femoral access only.

The recommended devices required for the retrieval of the OPTEASE Retrievable Filter are

listed in Table III below.

The IFU for the retrieval of the OPTEASE Retrievable Filter is contained in the packaging of the Cordis OPTEASE Retrieval Catheter.

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The OPTEASE Retrievable Filter can be retrieved up to and including 12 days after implantation from femoral access only. The recommended devices required for the retrieval of the OPTEASE Retrievable Filter are listed in Table III below. The IFU for the retrieval of the OPTEASE Retrievable Filter is contained in the packaging of the Cordis OPTEASE Retrieval Catheter.

Table III. Recommended Devices for the Retrieval of the OPTEASE® Retrivable Filter

Recommended Devices	Dimension	Cat. No.
OPTEASE® Retrieval Catheter with radiopaque tip or equivalent	10F, 80 cm length	466-C210F
Cordis J-tip EMERALD® 0.035" (0.89 mm) Guidewire or equivalent	260 cm length	502-455
Cordis AVANTI®+ Catheter Sheath Introducer or equivalent	10F (3.3 mm), 11 cm length	504-610X
Goose Neck Snare Kit or equivalent	10 mm loop diameter, 4F catheter	GN1000
	15 mm loop diameter, 6F catheter	GN1500
	20 mm loop diameter, 6F catheter	GN2000
	25 mm loop diameter, 6F catheter	GN2500
	30 mm loop diameter, 6F catheter	GN3000
All 120 cm snare length and 102 cm catheter length		

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- The OPTEASE Retrievable Filter can be retrieved up to and including 12 days after placement. The OPTEASE Retrievable Filter is considered a *permanent* implant if it is not retrieved within the specified time period.

- Retrieval is performed using the Cordis OPTEASE Retrieval Catheter or equivalent, and the recommended accessories (Table III). These devices used for filter retrieval are not included in the OPTEASE Retrievable Filter Introduction Set. The IFU for OPTEASE Retrievable Filter retrieval is contained in the Cordis OPTEASE Retrieval Catheter packaging.

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- The OPTEASE Retrievable Filter is supplied constrained in a plastic storage tube indicating the appropriate orientation for jugular approach. Never reload a fully ejected filter into the storage tube as this could result in incorrect orientation for the selected access site. Never reload a (partially) ejected filter into the storage tube as this could affect its shape, and function. Accordingly, Cordis will not be responsible for any direct, incidental or consequential damages resulting from replacement of the OPTEASE Retrievable Filter in the plastic storage tube.

- Filter fracture is a potential complication of vena cava filters. Anatomic locations that create concentrated stress points from filter deformation (for example, deployment at apex of scoliosis, overlapping of either of

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the renal ostia, or placement adjacent to a vertebral osteophyte) may contribute to fracture of a particular filter strut. However, reports of adverse clinical sequelae from filter fractures are rare.

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- Store in a cool, dark, dry place.
- Do not use if package is open or damaged.
- Use prior to "Use By" date.
- Do not autoclave or resterilize by any other method. Do not expose to temperatures above 54°C (130°F) for any length of time.
- This product is designed and intended for single use. It is not designed to undergo reprocessing and re-sterilization after initial use. Reuse of this product, including after reprocessing and/or re-sterilization, may cause a loss of structural integrity which could lead

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to a failure of the device to perform as intended and may lead to a loss of critical labeling/use information all of which present a potential risk to patient safety.

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- The safety and effectiveness of the OPTEASER Filter has not been established in patients who are known to be pregnant or who are lactating.

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- Normal operating mode of the MR system. The OPTEASE Retrievable Filter should not migrate in this MRI environment. Non-clinical testing has not been performed to rule out the possibility of filter migration at field strengths higher than 3 Tesla.

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Introducer by releasing the snap-fit ring at the hub (Figure 5). Withdraw the Angiographic Vessel Dilator.

Caution: To avoid damage to the Sheath Introducer tip, do not withdraw the dilator until the Sheath Introducer tip is at the desired location in the IVC.

Caution: Do not leave a Sheath Introducer in place for extended periods of time without a dilator or an obturator to support the cannula wall.

14. Aspirate from the sideport extension to remove any potential air.

15. According to the selected venous access site, determine which end of the storage tube (containing the filter) is to be placed into the Sheath Introducer hemostasis valve. This is indicated by the printed colored arrows and text (femoral: green; jugular/ antecubital: blue) on the storage tube. The arrow of the desired access site will point into the Sheath Introducer hemostasis valve (Figure 6A/B).

а. **NOTE:** The OPTEASE Filter is indicated for retrieval from the femoral approach only. The retrieval hook should be positioned caudally to allow the anti-migration fixation hooks to engage the caval wall properly.

16. Slowly advance the filter into the Sheath Introducer by advancing the obturator through the end of the storage tube (Figure 6B) until the filter is positioned well into the cannula of the Sheath Introducer.

17. To restore hemostasis, withdraw the storage tube out of the Sheath Introducer hemostasis valve. The storage tube can be fixed on the hub-end of the obturator by sliding it over the bump (Figure 7A).

18. Continue to advance the filter until the marker on the obturator is positioned at the Sheath Introducer hemostasis valve. This indicates that the filter is at the distal tip of the Sheath Introducer but still fully within the Sheath Introducer (Figure 7B).

19. Filter release (or deployment) is performed under continuous fluoroscopy. Prior to releasing the OPTEASE Retrievable Filter from the Sheath Introducer ensure that:

а. the intended filter location in the inferior vena cava is correct. Should the filter loca-

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tion be incorrect, reposition the Sheath Introducer; and b. the filter retrieval hook is oriented towards the caudal end of the vena cava. Should the filter orientation be incorrect, remove the Sheath Introducer (with filter inside) from the patient and discard the system. Recommence the procedure using a new sterile Sheath Introducer and storage tube with filter.

Note: The constrained, unexpanded, length of the OPTEASE Retrievable Filter is 67 mm. During release from the Sheath Introducer, the length of the filter will shorten as it expands. The ends of the filter will shorten equally as they move toward the axial center of the filter. Filter shortening is depicted in Figure 2. The length of the deployed filter is determined by its expanded diameter in the IVC. The shortening at each end of the expanded filter and the corresponding filter length is provided in Table II.

The position of the expanded filter within the IVC can be determined by selecting either the midpoint (i.e., axial center) or one end of the constrained filter as a reference point. When the filter is positioned by using the midpoint of the constrained filter as a reference point, the position of the midpoint of the expanded filter will not change. When the filter is positioned by selecting one end of the constrained filter as a reference point, the end of the expanded filter will move (i.e., filter shortening) toward the axial center.

20. To release the filter, keep the obturator fixed in position, then pull the Sheath Introducer

back over the obturator (Figure 8A).

21. Ensure that the filter is fully released and deployed (Figure 8B).

22. Once the OPTEASER Retrievable Filter has been deployed, remove the obturator from the Sheath Introducer.

23. Perform a control cavogram before terminating the procedure.

24. Remove the Sheath Introducer when clinically indicated by placing compression on the

vessel above the puncture site, and slowly withdrawing the Sheath Introducer.

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The OPTease Retrievable Filter can be retrieved up to and including 12 days after implantation from femoral access only.

The recommended devices required for the retrieval of the OPTease Retrievable Filter are listed in Table III below.

The IFU for the retrieval of the OPTease Retrievable Filter is contained in the packaging of the Cordis OPTease Retrieval Catheter.

Table III. Recommended Devices for the Retrieval of the OPTease® Retrivable Filter

Recommended Devices	Dimension	Cat. No.
OPTease® Retrieval		
Catheter with radiopaque tip or equivalent	10F, 80 cm length	466-C210F
Cordis J-tip EMERALD®		
0.035" (0.89 mm) Guidewire or equivalent	260 cm length	502-455
Cordis AVANTI®+ Catheter		
Sheath Introducer or equivalent	10F (3.3 mm), 11 cm length	504-610X
Goose Neck Snare Kit		
or equivalent	10 mm loop diameter, 4F catheter	GN1000
	15 mm loop diameter, 6F catheter	GN1500
	20 mm loop diameter, 6F catheter	GN2000
	25 mm loop diameter, 6F catheter	GN2500
	30 mm loop diameter, 6F catheter	GN3000
	All 120 cm snare length and 102 cm catheter length	

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- Do not use with Ethiodol or Lipiodol contrast media (Ethiodol and Lipiodol are trademarks of Guerbet S.A.)
- Do not use with Ethiodol** or Lipiodol contrast media (Ethiodol and Lipiodol are trademarks of Guerbet S.A), or other such contrast media, which incorporate the components of these agents.
- Do not use Ethiodol (Ethiodol is a trademark of Guerbet S.A.) or Lipiodol contrast medium.