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CERTIFICATE

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of FLYCAT ELECTRICAL CO.,LTD on the annexed DOCUMENT is genuine.

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INSTRUCTION FOR USE

Irrigator oral dental stationery model Revyline RL 100

Name:

Irrigator oral dental stationery model Revyline RL 100 in composition в составе:

1. main unit – 1 pc
2. water tank – 1 pc
3. standard nozzle – 3 pc
4. periodontal nozzle - 1 pc
5. orthodontic nozzle-1 pc
6. tongue cleaning nozzle-1 pc
7. implant nozzle -1 pc
8. operating instructions – 1 pc

(hereinafter-irrigator, product, device, , Revyline RL 100)

Thank you for purchasing the Revyline RL 100 oral irrigator. The use of an irrigator significantly increases the level of hygienic protection of the oral cavity, prevents the occurrence and development of dental diseases. The pulsating water flow stimulates the blood supply to the gums and improves their condition. The combination of high pressure and pulsating water supply technology helps the Revyline RL 100 irrigator to clean the interdental spaces, the cervical area of the teeth and the gingival furrow, which is not available for a regular toothbrush.

Intended use

Irrigators oral dental with accessories are designed to prevent diseases of the teeth and gums.

Indications

Irrigators oral dental with accessories indicated for:

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Precautions

Before using the Revyline RL 100 irrigator, read the operating instructions carefully and follow the safety instructions. This will help you use the device as efficiently and comfortably as possible.

1. Do not handle the plug with wet hands.
2. Do not allow the device body to come into contact with water or other liquids.
3. Do not use the irrigator while showering.
4. If the appliance falls into the water, do not touch it and immediately unplug it.
5. If the device was in a negative temperature environment, allow it to warm up to room temperature for 12 hours before use.
6. Always check the integrity of the power cable before each use of the appliance.
7. Connect the appliance only to an electrical outlet designed for the appropriate load.
8. Do not direct the water jet under the tongue, into the ear, nose cavity, or other vulnerable areas. The powerful water jet generated by the device can disrupt the integrity of the skin and mucous membranes in these areas.
9. Use the device in accordance with the instructions for use and/or recommendations of the dentist.
10. Use additional accessories and attachments only from the Revyline RL 100 irrigator.
11. Do not insert foreign objects into the connectors and holes on the device body.
12. Do not add granular substances, such as salt, to the water for the irrigator. They clog the pump and the handle of the irrigator.
13. Do not use mouthwashes containing iodine, oil, or bleaching agents.
14. Use water or a solution recommended by your dentist for irrigation.
15. Do not use the device if there are open wounds in the mouth or on the tongue.
16. If you are at risk for health reasons and are taking antibiotics before dental treatment, consult your dentist before using the irrigator.

17. Do not remove the rubber tip from the transparent nozzle for cleaning periodontal pockets. If it is lost, replace the nozzle with a new one.
18. Do not allow children to play with the device.

Operation of the device

1. Plug the power cord into a properly functioning outlet.
2. Remove the tank cover and the tank itself from the instrument body. Fill the tank with water. Secure the tank tightly to the instrument body.
3. Insert the nozzle into the central connector located at the top of the handle. The marking ring of the nozzle should be at the same level with the upper border of the handle.
4. To remove the nozzle, press and hold the nozzle lock button while removing it from the handle.
5. Bend down to the sink, insert the irrigator nozzle into the oral cavity and place its tip opposite the tooth. To turn on the device, turn the pressure power regulator clockwise. To select the water jet head power, use the pressure regulator: clockwise to increase, counterclockwise to decrease. Choose the most comfortable pressure level for yourself. During irrigation, move the nozzle along the gum line, consistently cleaning each tooth.
6. After using the device, turn the liquid pressure regulator counterclockwise until it clicks.

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Pause mode

The pause button located on the handle of the irrigator stops the water supply in the off position. It is forbidden to leave the device on for more than 10 seconds with the button in the off position, as this may cause it to break.

Built-in timer

The device is equipped with an automatic timer. Disconnection occurs after 3 minutes.

Selection the pressure level mode

The pressure level of the irrigation flow is set by a rotary regulator located on the device body. The design of the device allows you to choose one of ten intensity modes.

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
Attention! Some antiseptic solutions may contain insoluble particles, which, when in contact with the irrigation system, can affect its performance and reduce its service life. After each use of the special solution for the irrigator, rinse the system thoroughly. Fill the tank with enough water, point the nozzle in the direction of the sink, and turn on the device. Run all the water in the tank through the irrigation system. At the same time, the remaining solutions are washed out of the irrigator, protecting the system from clogging.

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Cleaning and maintenance

Before cleaning the irrigator, disconnect the power supply from the electrical outlet. If necessary, use a soft cloth and a mild, non-abrasive detergent for cleaning. Do not immerse the nozzles in boiling water. To clean them, use alcohol. With regular use of hard water, calcareous deposits may form on the parts of the device. To eliminate them, use a special Revelyne product for cleaning the irrigator, according to the instructions for its use. The manufacturer recommends cleaning the irrigator from one to three times every three months, depending on the hardness of the water.

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Operational temperature	from +15°C to +30°C
Storage temperature	from +15°C to +30°C
Relative humidity	from 10% to 90%.
Atmospheric pressure	N/A

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Volume	600 ml
Pressure modes	10
Ripple frequency	1200 – 1700 per min
Housing moisture protection	IPX7
Supply	AC power supply 100-240V, 50/60 Hz
Operating time at full tank	90 sec
Maximum continuous use time (Auto-Off)	3 min
Power consumption, no more than	18 Wt
Noise level	≤ 78 dB

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Solution: After turning on the device, use the pressure regulator to select the desired level of jet power, which is selected individually. To increase the efficiency of irrigation, use a water pressure that does not cause you discomfort.

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic emission test	Compliance	Electromagnetic environment-instructions
Radio interference by CISPR 11	Group 1	Dental oral irrigators with accessories use radio frequency energy only to perform internal functions. The level of radio frequency interference emission is low and is not likely to cause disruption to nearby electronic equipment.
Radio interference by CISPR 11	Class B	
Harmonic components of the current according to IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Comply	

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.			
Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Electrostatic discharges (ESR) according to IEC 61000-4-2	Contact discharge ± 6 kV	Contact discharge ± 6 kV	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
	Air discharge ± 8 kV	Air discharge ± 8 kV	
Nanosecond pulse interference according to IEC 61000-4-4	For ± 2 kV power supply lines For input/ output lines ± 1 kV	For power supply lines of ± 2 kV	The quality of electrical energy in the network in accordance with the typical conditions of a commercial or hospital environment.
Microsecond high-energy pulse interference according to IEC 61000-4-5	When applying interference according to the "wire-to-wire" scheme ± 1 kV When applying interference according to the	When applying interference according to the "wire-to-wire" scheme ± 1 kV	The quality of the power supply should be consistent with a standard commercial or hospital environment.

	"wire-to-ground" scheme ± 2 kV		
voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	The quality of electrical energy in the network is in accordance with the typical conditions of a commercial or hospital environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.

NOTE: U_H – the voltage level of the electrical network until the moment of applying the test action.

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.			
Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Conductive interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 Radio frequency electromagnetic field according to IEC 61000-4-3	3 V (mean-square value) in the band from 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any element of the Irrigator, including cables, should not be less than the recommended spatial separation, which is calculated in accordance with the following expressions for the frequency of the transmitter. Recommended spatial spacing: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ (80 MHz to 800 MHz) $d=2.3\sqrt{P}$ (800 MHz to 2.5 GHz), where d is the recommended spatial spacing, m; P – the rated maximum output power of the transmitter, the maximum output power of the transmitter, W, set by the manufacturer. The field strength in the distribution of radio waves from stationary radio transmitters, according to the results of observations of the electromagnetic environment a) should be lower than the level of compliance in each of the frequencies b). Interference may occur in the vicinity of equipment marked with the sign: 

NOTE 1: For frequencies between 80 MHz and 800 MHz, a higher field strength value is applied.
NOTE 2: Expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people.
a: The field strength during the propagation of radio waves from stationary radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and terrestrial stationary radio stations, amateur radio stations, AM and FM radio transmitters, television transmitters cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the Irrigator exceed the applicable compliance levels, the operation of the Irrigator should be monitored to verify their normal functioning. If the monitoring process reveals a deviation from normal operation, then additional measures may need to be taken, such as reorienting or moving the Irrigator. b: Outside the 150 kHz to 80 MHz band, the field strength must be less than 3 V / m

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Irrigators are intended for use in an electromagnetic environment in which the levels of radiated interference are monitored. The buyer or user of irrigators can avoid the influence of electromagnetic interference by ensuring a minimum spatial separation between portable and mobile radio frequency communication devices (transmitters) and irrigators as recommended below, taking into account the maximum output power of the communication devices.			
Rated maximum output power of the transmitter P, W	Spatial spacing d, m, depending on the frequency of the transmitter		
	d = $1.2\sqrt{P}$ in the band from 150 kHz to 80 MHz	d = $1.2\sqrt{P}$ in the band from 80 to 800 MHz	d = $2.3\sqrt{P}$ in the band from 800 MHz to 1.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23
When determining the recommended values of spatial spacing d for transmitters with a nominal maximum output power not specified in the table, the given expressions are substituted for the nominal maximum output power P in watts, specified in the documentation of the transmitter manufacturer.			
NOTE 1: At frequencies of 80 and 800 MHz, a higher field strength is applied.			
NOTE 2: These expressions do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people			

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AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE RUSSIAN FEDERATION

Obushenkova individual entrepreneur Maxim Alekseevich
Udmurtskaya republic, Izhevsk, Karl Marks street, bl 425, bl. 27, flat 6.
Tel. 8 (800) 775-38-49
www.revyline.ru.

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of FLYCAT ELECTRICAL CO., LTD on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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证明 网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Irrigator oral dental stationery model Revyline RL 500

Name:

Irrigator oral dental stationery mode Revyline RL 500 in composition:

1. main unit – 1 pc
2. water tank – 1 pc
3. standard nozzle – 3 pc
4. periodontal nozzle - 1 pc
5. orthodontic nozzle-1 pc
6. tongue cleaning nozzle-1 pc
7. implant nozzle -1 pc
8. instruction for use – 1 pc

Accessories:

1. Nozzle case- 1 pc

(hereinafter-irrigator, product, device, , Revyline RL 500)

Thank you for purchasing the Revyline RL 500 oral irrigator. The use of an irrigator significantly increases the level of hygienic protection of the oral cavity, prevents the occurrence and development of dental diseases. The pulsating water flow stimulates the blood supply to the gums and improves their condition. The combination of high pressure and pulsating water supply technology helps the Revyline RL 500 irrigator to clean the interdental spaces, the cervical area of the teeth and the gingival furrow, which is not available for a regular toothbrush.

Intended use

Irrigators oral dental are designed to prevent diseases of the teeth and gums.

Indications

Irrigators oral dental indicated for:

- prevention of bleeding gums, caries, periodontal diseases, gingivitis, peri-implantitis;
- comprehensive care and hygiene of the oral cavity and interdental spaces;
- removal of plaque, prevention of its mineralization and formation of tartar;
- oral care in the presence of bridges, implants, crowns and braces.

Contraindications

Currently, no contraindications have been identified

Precautions

Before using the Revyline RL 500 irrigator, read the operating instructions carefully and follow the safety instructions. This will help you use the device as efficiently and comfortably as possible.

1. Do not handle the plug with wet hands.
2. Do not allow the device body to come into contact with water or other liquids.
3. Do not use the irrigator while showering.
4. If the appliance falls into the water, do not touch it and immediately unplug it.
5. If the device was in a negative temperature environment, a room temperature for 12 hours before use.
6. Always check the integrity of the power cable before each use.
7. Connect the appliance only to an electrical outlet designed for the device.
8. Do not direct the water jet under the tongue, into the ear, nose cavity, or other vulnerable areas. The powerful water jet generated by the device can disrupt the integrity of the skin and mucous membranes in these areas.
9. Use the device in accordance with the instructions for use and/or recommendations of the dentist.
10. Use additional accessories and attachments only from the Revyline RL 100 irrigator.
11. Do not insert foreign objects into the connectors and holes on the device body.
12. Do not add granular substances, such as salt, to the water for the irrigator. They clog the pump and the handle of the irrigator.
13. Do not use mouthwashes containing iodine, oil, or bleaching agents.

- Use water or a solution recommended by your dentist for irrigation.
15. Do not use the device if there are open wounds in the mouth or on the tongue.
 16. If you are at risk for health reasons and are taking antibiotics before dental treatment, consult your dentist before using the irrigator.
 17. Do not remove the rubber tip from the transparent nozzle for cleaning periodontal pockets. If it is lost, replace the nozzle with a new one.
 18. Do not allow children to play with the device.

Operation of the device

1. Plug the power cord into a properly functioning outlet.
2. Remove the tank cover and the tank itself from the instrument body. Fill the tank with water. Secure the tank tightly to the instrument body.
3. Insert the nozzle into the central connector located at the top of the handle. The marking ring of the nozzle should be at the same level with the upper border of the handle.
4. To remove the nozzle, press and hold the nozzle lock button while removing it from the handle.
5. To turn on the device, press the power on/off button located on the pressure power regulator. To select the water jet head power, use the pressure adjustment knob located on the device body. Bend down to the sink, insert the irrigator nozzle into the oral cavity and place its tip opposite the tooth. During irrigation, move the nozzle along the gum line, consistently cleaning each tooth.
6. After using the device, turn the liquid pressure regulator to the minimum level, and then press the power on/off button.

Irrigation mode

Position the nozzle in the oral cavity so that the flow of the jet is directed at an angle of 45° to the surface of the tooth. To prevent water from splashing, close your lips loosely. During irrigation, water should flow naturally from the mouth into the sink. For best results, first clean the side teeth (molars), and then the front teeth. Move the tip of the nozzle along the gum line, holding it between the teeth. Continue the procedure until you clean the teeth inside and out, as well as all the interdental spaces.

Pause mode

The operation of the device can be stopped at any time by shifting the button located on the handle of the irrigator. The reverse shift of the button resumes the operation of the device.

Attention! The device should not be used for a long time in this mode.

Massage mode

To activate the mode, press the Massage button located above the pressure level control. When performing a massage, it is recommended to use a standard replaceable nozzle.

Selecting the pressure level mode

The pressure level of the irrigation flow is set by a rotary regulator located on the device body. The design of the device allows you to choose one of the seventeen intensity modes.

Attention! At the beginning of the application of the device, it is recommended to use a minimum pressure level. This will help you adapt to irrigation cleaning and avoid the risk of gum inflammation and bleeding.

After irrigation

Turn off the irrigator. Unplug the power cord from the wall outlet. The irrigator is not intended for storing water in it. After use, the remaining liquid must be drained.

UV disinfection of replaceable nozzles

The irrigator is equipped with the function of disinfection of the nozzles by ultraviolet irradiation. The lamp is installed in the nozzle storage compartment. To perform disinfection, the nozzles must be placed in the compartment and the compartment door must be tightly closed. After 15 minutes, the lamp turns off automatically. During the operation of the lamp under the influence of ultraviolet rays, atomic oxygen (ozone) is formed in the compartment, which destroys 99% of bacteria and viruses.

Attention! When the UV lamp is running, be careful not to get radiation in your eyes. High-frequency light waves can cause inflammation of the retina.

Use of mouthwashes and other solutions

For the best effect, it is recommended to use **Revyline** rinses and specialized balms. Calculate the dosage according to the instructions of the product used. After each use of the solutions, it is necessary to rinse the irrigator.

Some antiseptic solutions may contain insoluble particles, which, when in contact with the irrigation system, can affect its performance and reduce its service life. After each use of the special solution for the irrigator, rinse the system thoroughly. Fill the tank with enough water, point the nozzle in the direction of the sink, and turn on the device. Run all the water in the tank through the irrigation system. At the same time, the remaining solutions are washed out of the irrigator, protecting the system from clogging.

List of applicable standards

UL 1431	Personal hygiene products
ISO 9001:2015	Quality management systems. Requirements

Symbols used on labeling

Symbol	Description
	Manufacturer
	Date of manufacture (month, year)
	Caution! Special utilization required!
	Working part type BF
	Device class II
	Sign of product circulation on the territory of markets of participation of Custom union (in accordance with TR TS 020/2011)

Cleaning and maintenance

Before cleaning the irrigator, disconnect the power supply from the electrical outlet. If necessary, use a soft cloth and a mild, non-abrasive detergent for cleaning. Do not immerse the nozzles in boiling water. To clean them, use alcohol. With regular use of hard water, calcareous deposits may form on the parts of the device. To eliminate them, use a special Revyline product for cleaning the irrigator, according to the instructions for its use. The manufacturer recommends cleaning the irrigator from one to three times every three months, depending on the hardness of the water.

Transportation conditions

Irrigators oral dental transported by all types of transport in accordance with the rules of transportation applicable to this type of transport

Temperature	Relative humidity	Atmospheric pressure
from -20°C to +50°C	from 10% to 90%.	N/A

Operational and storage conditions

Operational temperature	from +15°C to +30°C
Storage temperature	from +15°C to +30°C

Relative humidity	from 10% to 90%.
Atmospheric pressure	N/A

Service life

The average service life of irrigators is 5 years.

Guaranteed by the manufacturer values of the main parameters, characteristics (properties) of the medical device/technical characteristics of the device

Model	Revyline RL 500
Parameter*	
Уровень давления жидкости	from 210 to 870 kPa
Tank volume	600 ml
Jet pressure modes	10
Ripple frequency	1200 – 1700 per min
Housing moisture protection	IPX7
Supply	AC power supply 100-240V, 50/60 Hz
Operating time at full tank	80 sec
Maximum continuous use time (Auto-Off)	3 min
Power consumption, no more than	18 Vt
Noise level	≤ 78 dB

* where applicable, the deviation from the specified characteristics is +/- 10%, unless otherwise specified

Requirements for the maintenance and repair of the medical device

Warranty maintenance of the device is carried out in the service center during the current warranty period (2 years free of charge). Further maintenance of the irrigator is carried out at the expense of the buyer.

Disposal of a medical device

The products should be disposed of in accordance with accepted regulations and should not be disposed of together with household waste.

Warranty obligations

The device is guaranteed for 2 years. Maintenance of the irrigator is carried out in the service center during the current warranty period free of charge. Further, the service of the device is performed at the expense of the buyer.

The warranty does not apply:

For consumables and accessories;

For goods that have failed and / or received defects due to their improper use; careless and / or negligent use (including overloading), resulting in damage; modification, opening and / or repair by an unauthorized organization (private person);

For products that have failed and / or become defective due to the connection to the wrong power source;

For goods that were operated with unremunerated defects;

For malfunctions that are caused by reasons beyond the control of the manufacturer, such as the actions of third parties, natural phenomena and natural disasters, domestic and wild animals, insects, ingestion of foreign objects and liquids;

For external and internal dirt, scratches, cracks, scuffs and other mechanical damage that occurred during operation and as a result of natural wear and tear.

Common problems and solutions

The nozzle is not removed.

Solution: While removing the nozzle, press and hold the nozzle lock button located on the instrument body

Difficulties in selecting the jet pressure level.

on: After turning on the device, use the pressure regulator to select the desired level of jet power, which is selected individually. To increase the efficiency of irrigation, use a water pressure that does not cause you discomfort.

Bleeding during irrigation.

Solution: Bleeding is a characteristic symptom of gum inflammation. First of all, it is necessary to reduce the pressure level of the jet. The irrigator removes plaque and bacteria from the gums. With further use, the bleeding will stop. For inflamed and weak gums, use a periodontal nozzle.

The maximum water pressure is uncomfortable.

Solution: The average liquid pressure level is suitable for most people. A comfortable liquid pressure can be adjusted by adjusting the liquid pressure level.

How to prevent water splashing?

Solution: During irrigation, do not stand upright with your mouth wide open. Lean towards the sink, open your mouth slightly. The water will naturally flow out of the mouth.

The engine is running, but no water is supplied.

Possible solutions:

1. Remove the water tank and re-attach it to the appliance body;
2. Increase the pressure level with the liquid pressure regulator;
3. The "Pause" switch can be in the OFF position. Check the status of the switch.

Leakage between the nozzle and the housing.

Solution: Check whether the nozzle is fully inserted. If necessary, remove the nozzle and then insert it again.

Features

The irrigator is equipped with a lid-pencil case, in which you can store the nozzles.

For reliable fixation of the irrigator to the horizontal surface, there are suction cups on the bottom of the device.

Information about electromagnetic compatibility

MANUFACTURER'S Manual and Declaration-ELECTROMAGNETIC EMISSION

Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic emission test	Compliance	Electromagnetic environment-instructions
Radio interference by CISPR 11	Group 1	Dental oral irrigators with accessories use radio frequency energy only to perform internal functions. The level of radio frequency interference emission is low and is not likely to cause disruption to nearby electronic equipment.
Radio interference by CISPR 11	Class B	
Harmonic components of the current according to IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Comply	Dental oral irrigators with accessories are suitable for use in all institutions, including residential premises and premises directly connected to the low-voltage power supply network that provides buildings used for domestic purposes.

MANUFACTURER'S Manual and Declaration-NOISE IMMUNITY

Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.			
Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Electrostatic discharges (ESR) according to IEC 61000-4-2	Contact discharge ± 6 kV	Contact discharge ± 6 kV	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the
	Air discharge ± 8 kV	Air discharge ± 8 kV	

			relative humidity should be at least 30%.
nanosecond pulse interference according to IEC 61000-4-4	For ± 2 kV power supply lines For input/ output lines ± 1 kV	For power supply lines of ± 2 kV	The quality of electrical energy in the network in accordance with the typical conditions of a commercial or hospital environment .
Microsecond high-energy pulse interference according to IEC 61000-4-5	When applying interference according to the "wire-to-wire" scheme ± 1 kV When applying interference according to the "wire-to-ground" scheme ± 2 kV	When applying interference according to the "wire-to-wire" scheme ± 1 kV	The quality of the power supply should be consistent with a standard commercial or hospital environment.
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11	<5% U_H (voltage drop >95% U_H) for 0.5 period. 40% U_H (voltage drop >60% U_H) for 5 periods. 70% U_H (voltage drop >30% U_H) for 25 periods. <5% U_H (>95% voltage drop U_H) for 5 seconds	<5% U_H (voltage drop >95% U_H) for 0.5 period. 40% U_H (voltage drop >60% U_H) for 5 periods. 70% U_H (voltage drop >30% U_H) for 25 periods. <5% U_H (>95% voltage drop U_H) for 5 seconds	The quality of electrical energy in the network is in accordance with the typical conditions of a commercial or hospital environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.
NOTE: U_H – the voltage level of the electrical network until the moment of applying the test action.			

MANUFACTURER'S Manual and Declaration-NOISE IMMUNITY-for non-LIFE-SUPPORT medical devices

Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.

Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Conductive interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6	3 V (mean-square value) in the band from 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any element of the Irrigator, including cables, should not be less than the recommended spatial separation, which is calculated in accordance with the following expressions for the frequency of the transmitter. Recommended spatial spacing: $d=1.2\sqrt{P}$

Radio frequency electromagnetic field according to IEC 61000-4-3			$d=1.2\sqrt{P}$ (80 MHz to 800 MHz) $d=2.3\sqrt{P}$ (800 MHz to 2.5 GHz), where d is the recommended spatial spacing, m; P – the rated maximum output power of the transmitter, the maximum output power of the transmitter, W, set by the manufacturer. The field strength in the distribution of radio waves from stationary radio transmitters, according to the results of observations of the electromagnetic environment a) should be lower than the level of compliance in each of the frequencies b). Interference may occur in the vicinity of equipment marked with the sign: 
--	--	--	--

Note 1: For frequencies between 80 MHz and 800 MHz, a higher field strength value is applied.

Note 2: Expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people.

a: The field strength during the propagation of radio waves from stationary radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and terrestrial stationary radio stations, amateur radio stations, AM and FM radio transmitters, television transmitters cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the Irrigator exceed the applicable compliance levels, the operation of the Irrigator should be monitored to verify their normal functioning. If the monitoring process reveals a deviation from normal operation, then additional measures may need to be taken, such as reorienting or moving the Irrigator. b: Outside the 150 kHz to 80 MHz band, the field strength must be less than 3 V / m

Recommended values of spatial separation between portable and mobile radio frequency communication devices and oral irrigators dental with accessories

Irrigators are intended for use in an electromagnetic environment in which the levels of radiated interference are monitored. The buyer or user of irrigators can avoid the influence of electromagnetic interference by ensuring a minimum spatial separation between portable and mobile radio frequency communication devices (transmitters) and irrigators as recommended below, taking into account the maximum output power of the communication devices.

Rated maximum output power of the transmitter P, W	Spatial spacing d, m, depending on the frequency of the transmitter		
	d = $1.2\sqrt{P}$ in the band from 150 kHz to 80 MHz	d = $1.2\sqrt{P}$ in the band from 80 to 800 MHz	d = $2.3\sqrt{P}$ in the band from 800 MHz to 1.5 GHz
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

When determining the recommended values of spatial spacing d for transmitters with a nominal maximum output power not specified in the table, the given expressions are substituted for the nominal maximum output power P in watts, specified in the documentation of the transmitter manufacturer

NOTE 1: At frequencies of 80 and 800 MHz, a higher field strength is applied.

NOTE 2: These expressions do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people

Manufacturer

Flycat Electrical Co., Ltd., 101-401 of No.28 Bldg & 101, 301, 401 of 30 Bldg, Longwangmiao Industrial Park, Baishixia Community, Fuyong Subdistrict, Bao'an District, Shenzhen City, Guangdong Province, China

**AUTHORIZED REPRESENTATIVE OF THE MANUFACTURER IN THE TERRITORY OF THE
RUSSIAN FEDERATION**

Obushenkova individual entrepreneur Maxim Alekseevich
Udmurtskaya republic, Izhevsk, Karl Marks street, bl. 425, bl. 27, flat 6.
Tel. 8 (800) 775-38-49
www.revyline.ru.

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证 询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Irrigator oral dental stationery model Revyline RL 900

Name:

Irrigator oral dental stationery model Revyline RL 900 in composition в составе:

1. main unit – 1 pc
2. water tank – 1 pc
3. standard nozzle – 3 pc
4. periodontal nozzle - 1 pc
5. orthodontic nozzle-1 pc
6. tongue cleaning nozzle-1 pc
7. implant nozzle -1 pc
8. instruction for use – 1 pc

(hereinafter-irrigator, product, device, , Revyline RL 900)

Dear customer, you have become the owner of a high-tech device that will allow your teeth and gums to be beautiful and healthy. Experience a new level of purity after the first application.

Important!!!

Before using the Revyline RL 900 irrigator, read the operating instructions carefully and follow the safety rules. This will help you use the device as efficiently and comfortably as possible.

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Precautions

Before using the Revyline irrigator, read the operating instructions carefully and follow the safety instructions. This will help you use the device as efficiently and comfortably as possible.

1. Do not handle the plug with wet hands.
2. Do not allow the device body to come into contact with water or other liquids.
3. Do not use the irrigator while showering.
4. If the appliance falls into the water, do not touch it and immediately pull the plug out of the socket.
5. If the device was in a negative temperature environment, allow it to warm up to room temperature for 12 hours before use.
6. Always check the integrity of the power cable before each use of the appliance.
7. Connect the appliance only to an electrical outlet designed for the appropriate load.
8. Do not direct the water flow under the tongue, into the ears, nose, or other sensitive areas. The device generates a powerful pulsating water wave, which can disrupt the integrity of the skin and mucous membranes.
9. Use the device in accordance with the instructions for use and/or recommendations of the dentist.
10. Use only branded attachments and accessories.
11. Do not insert foreign objects into the connectors and holes on the device body.
12. Do not add granular substances, such as salt, to the water for the irrigator. They clog the pump and the handle of the irrigator.
13. Do not fill the tank with rinses containing iodine, oil, bleaching ingredients, as well as rinses that allow natural precipitation and rinses that do not completely dissolve in water. They reduce the efficiency and service life of the device.
14. Use water or a solution recommended by your dentist for irrigation.
15. Do not use the device if there are open wounds in the mouth or on the tongue.

- If you are at risk for health reasons and take antibiotics before performing dental treatment, consult your dentist before using the irrigator.
17. Do not remove the rubber tip from the transparent nozzle for cleaning periodontal pockets. If it is lost, replace the nozzle with a new one.
 18. Do not allow children to play with the appliance.
 19. The use of the irrigator by children is possible from the age of 6 only under the supervision of adults.

Operation of the device

1. Plug the power cord into a properly functioning outlet.
2. Remove the case-cover, remove the tank from the case, fill it with warm purified water. Install the tank inside the instrument housing.
3. Insert the nozzle into the connector located on the top of the handle and insert it until it clicks into place. To remove the nozzle, hold down the lock button and remove the nozzle from the device body. Attention! Do not press the eject button while the device is in operation.
4. To turn on the device, press the power on/off button located on the handle. To select the water jet pressure, use the pressure control buttons (+/-) located on the handle of the device below the power button. Bend down to the sink, insert the irrigator nozzle into the oral cavity and place its tip opposite the tooth. During irrigation, move the nozzle along the gum line, consistently cleaning each tooth.
5. After using the device, turn the liquid pressure regulator to the minimum level, and then press the power on/off button.

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The selected pressure mode is displayed on the electronic display.

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
	Manufacturer
	Date of manufacture (month, year)
	Caution! Special utilization required!
	Working part type BF
	Device class II
	Sign of product circulation on the territory of markets of participation of Custom union (in accordance with TR TS 020/2011)

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ply	AC power supply 100-240V, 50/60 Hz
Operating time at full tank	76 sec
Maximum continuous use time (Auto-Off)	3 min
Power consumption, no more than	18 Vt
Noise level	≤ 78 dB

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Possible solutions:

1. Remove the water tank and re-attach it to the appliance body;
2. Increase the pressure level with the liquid pressure regulator;

The "Pause" switch can be in the OFF position. Check the status of the switch.

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic emission test	Compliance	Electromagnetic environment-instructions
Radio interference by CISPR 11	Group 1	Dental oral irrigators with accessories use radio frequency energy only to perform internal functions. The level of radio frequency interference emission is low and is not likely to cause disruption to nearby electronic equipment. Dental oral irrigators with accessories are suitable for use in all institutions, including residential premises and premises directly connected to the low-voltage power supply network that provides buildings used for domestic purposes.
Radio interference by CISPR 11	Class B	
Harmonic components of the current according to IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Comply	

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.			
Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Electrostatic discharges (ESR) according to IEC 61000-4-2	Contact discharge ± 6 kV	Contact discharge ± 6 kV	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
	Air discharge ± 8 kV	Air discharge ± 8 kV	
Nanosecond pulse interference according to IEC 61000-4-4	For ± 2 kV power supply lines For input/ output lines ± 1 kV	For power supply lines of ± 2 kV	The quality of electrical energy in the network in accordance with the typical conditions of a commercial or hospital environment.
Microsecond high-energy pulse interference according to IEC 61000-4-5	When applying interference according to the "wire-to-wire" scheme ± 1 kV When applying interference according to the "wire-to-ground" scheme ± 2 kV	When applying interference according to the "wire-to-wire" scheme ± 1 kV	The quality of the power supply should be consistent with a standard commercial or hospital environment.
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period.	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period.	The quality of electrical energy in the network is in accordance with the typical conditions of a commercial or hospital
	$40\% U_H$ (voltage	$40\% U_H$ (voltage	

	drop >60% U_H) for 5 periods. 70% U_H (voltage drop >30% U_H) for 25 periods. <5% U_H (>95% voltage drop U_H) for 5 seconds	drop >60% U_H) for 5 periods. 70% U_H (voltage drop >30% U_H) for 25 periods. <5% U_H (>95% voltage drop U_H) for 5 seconds	environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.
NOTE: U_H – the voltage level of the electrical network until the moment of applying the test action.			

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Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Conductive interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 Radio frequency electromagnetic field according to IEC 61000-4-3	3 V (mean-square value) in the band from 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any element of the Irrigator, including cables, should not be less than the recommended spatial separation, which is calculated in accordance with the following expressions for the frequency of the transmitter. Recommended spatial spacing: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ (80 MHz to 800 MHz) $d=2.3\sqrt{P}$ (800 MHz to 2.5 GHz), where d is the recommended spatial spacing, m; P – the rated maximum output power of the transmitter, the maximum output power of the transmitter, W, set by the manufacturer. The field strength in the distribution of radio waves from stationary radio transmitters, according to the results of observations of the electromagnetic environment a) should be lower than the level of compliance in each of the frequencies b). Interference may occur in the vicinity of equipment marked with the sign: 

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a: The field strength during the propagation of radio waves from stationary radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and terrestrial stationary radio stations, amateur radio stations, AM and FM radio transmitters, television transmitters cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the

Measured values at the location of the Irrigator exceed the applicable compliance levels, the operation of the Irrigator should be monitored to verify their normal functioning. If the monitoring process reveals a deviation from normal operation, then additional measures may need to be taken, such as reorienting or moving the Irrigator. b: Outside the 150 kHz to 80 MHz band, the field strength must be less than 3 V / m

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Obushenkova individual entrepreneur Maxim Alekseevich

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INSTRUCTION FOR USE

Irrigator oral dental portable model Revyline RL 210

Name:

Irrigator oral dental portable model Revyline RL 210 in composition:

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Accessories:

1. case – 1 pc
2. nozzle case – 1 pc
3. USB adapter – 1 pc
4. USB cable – 1 pc

(hereinafter-irrigator, product, device, Revyline RL 210)

Thank you for purchasing the Revyline RL 210 oral irrigator. The use of an irrigator significantly increases the level of hygienic protection of the oral cavity, prevents the occurrence and development of dental diseases. The pulsating water flow stimulates the blood supply to the gums and improves their condition. The combination of high pressure and pulsating water supply technology helps the Revyline RL 210 irrigator to clean the interdental spaces, the cervical area of the teeth and the gingival furrow, which is not available for a regular toothbrush.

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Precautions

Before using the Revyline RL 210 irrigator for the first time, read the operating instructions carefully and save them.

To prevent electrical injuries, burns, and fires:

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Turning on the device

Open the tank cover, tilt the body so that the inlet is horizontal and perpendicular to the direction of the jet. Fill the tank with boiled, purified water. Close the tank cover. To turn on the device, press the power ON/OFF button once. To select the desired operating mode of the device, use the MODE button. To turn off the irrigator, press the ON/OFF button again.

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Cleaning and maintenance

Do not immerse the device in water. It is forbidden to store liquid in the tank. Use a soft cloth and a soft, non-abrasive cleaner to clean the surface of the appliance. To facilitate the cleaning of the tank, it must be disconnected from the housing. With regular use of hard water, calcareous deposits may form on the parts of the device. To eliminate them, use a special Revelyne product for cleaning the irrigator, according to the instructions for its use. The manufacturer recommends cleaning the irrigator from one to three times every three months, depending on the hardness of the water.

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Model	Revyline RL 210
Parameter*	
Liquid pressure level	From 210 to 700 kPa
Tank volume	150 ml
Working modes	3
Ripple frequency	1400 – 1800 per minute
Housing moisture protection	IPX7
Supply	Lithium Battery 1400 mAh
Time of full charge	4 hours
Battery life	Till 500 cycles (charge/discharge)
Battery life on a single charge	Minimum 60 minutes
Operating time at full tank	At normal mode (Normal) 50 sec At soft mode (Soft) 75 sec At pulse massage mode (Pulse) 65 sec
Auto-Off	120 sec
Power consumption, no more than	5 Vt
Noise level	≤ 78 dB

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic emission test	Compliance	Electromagnetic environment-instructions
Radio interference by CISPR 11	Group 1	Dental oral irrigators with accessories use radio frequency energy only to perform internal functions. The level of radio frequency interference emission is low and is not likely to cause disruption to nearby electronic equipment.
Radio interference by CISPR 11	Class B	Dental oral irrigators with accessories are suitable for use in all institutions, including residential premises and premises directly connected to the low-voltage power supply network that provides buildings used for domestic purposes.
Harmonic components of the current according to IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Comply	

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Dental oral irrigators with accessories are designed for use in the electromagnetic environment defined below. The buyer or user of the product should ensure that it is used in the specified electromagnetic environment.			
Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Electrostatic discharges (ESR) according to IEC 61000-4-2	Contact discharge ± 6 kV Air discharge ± 8 kV	Contact discharge ± 6 kV Air discharge ± 8 kV	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
Nanosecond pulse interference according to IEC 61000-4-4	For ± 2 kV power supply lines For input/ output lines ± 1 kV	For power supply lines of ± 2 kV	The quality of electrical energy in the network in accordance with the typical conditions of a commercial or hospital environment.
Microsecond high-energy pulse interference according to IEC 61000-4-5	When applying interference according to the "wire-to-wire" scheme ± 1 kV When applying interference according to the "wire-to-ground" scheme ± 2 kV	When applying interference according to the "wire-to-wire" scheme ± 1 kV	The quality of the power supply should be consistent with a standard commercial or hospital environment.
Voltage dips, short-term interruptions and	$<5\%$ U_H (voltage drop $>95\%$ U_H) for	$<5\%$ U_H (voltage drop $>95\%$ U_H) for 0.5	The quality of electrical energy in the network is in

voltage changes in the input power lines according to IEC 61000-4-11	0.5 period. 40% U_H (voltage drop $>60\% U_H$) for 5 periods. 70% U_H (voltage drop $>30\% U_H$) for 25 periods. <5% U_H ($>95\%$ voltage drop U_H) for 5 seconds	period. 40% U_H (voltage drop $>60\% U_H$) for 5 periods. 70% U_H (voltage drop $>30\% U_H$) for 25 periods. <5% U_H ($>95\%$ voltage drop U_H) for 5 seconds	accordance with the typical conditions of a commercial or hospital environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.

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a: The field strength during the propagation of radio waves from stationary radio transmitters, such as

base stations of radiotelephone networks (cellular/wireless), and terrestrial stationary radio stations, amateur radio stations, AM and FM radio transmitters, television transmitters cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the Irrigator exceed the applicable compliance levels, the operation of the Irrigator should be monitored to verify their normal functioning. If the monitoring process reveals a deviation from normal operation, then additional measures may need to be taken, such as reorienting or moving the Irrigator. b: Outside the 150 kHz to 80 MHz band, the field strength must be less than 3 V / m

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证明一网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Irrigator oral dental portable model Revyline RL 200, Revyline RL 200 XL

Name:

Irrigator oral dental portable model Revyline RL 200, Revyline RL 200 XL in composition:

1. main unit – 1 pc
2. water tank – 1 pc
3. standard nozzle – 2 pc
4. instruction for use – 1 pc

Accessories:

1. case – 1 pc
2. nozzle case – 1
3. USB adapter – 1 pc
4. USB cable – 1 pc

(hereinafter-irrigator, product, device, Revyline RL 200, Revyline RL 200 XL)

Thank you for purchasing the Revyline RL 200 (Revyline RL 200 XL) oral irrigator. The use of an irrigator significantly increases the level of hygienic protection of the oral cavity, prevents the occurrence and development of dental diseases. The pulsating water flow stimulates the blood supply to the gums and improves their condition. The combination of high pressure and pulsating water supply technology helps the Revyline RL 200 (Revyline RL 200 XL) irrigator to clean the interdental spaces, the cervical area of the teeth and the gingival furrow, which is not available for a regular toothbrush.

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Precautions

Before using the Revyline RL 200 (Revyline RL 200 XL) irrigator, read the operating instructions carefully and follow the safety instructions. This will help you use the device as efficiently and comfortably as possible.

1. Do not immerse the appliance in water or other liquid.
2. Do not use the irrigator while showering.
3. Do not touch the charger connected to the mains with wet hands.
4. Do not store the appliance in places where it may fall into a water container or sink.
5. Do not use the device if the charger is faulty or the power cable is damaged.
6. Do not direct the water jet under the tongue, in the ears, nose cavity or other vulnerable areas. The powerful water jet generated by the device can disrupt the integrity of the skin and mucous membranes in these areas.
7. Use only branded attachments and accessories.
8. Do not insert foreign objects into the replacement nozzle socket
9. Charge the device away from hot surfaces.
10. Do not turn on the irrigator in a room where oxygen or aerosols have previously been used.
11. Do not fill the tank with rinses containing iodine, oil, or bleaching ingredients. They reduce the efficiency and service life of the device.
12. Fill the tank with water or a solution recommended by your dentist.
13. Before using the irrigator, remove all removable structures from the oral cavity.
14. Do not use the device if there are open wounds in the mouth or on the tongue.
15. If you are at risk for health reasons and are taking antibiotics before dental treatment, consult your dentist before using the irrigator.

16. The use of the irrigator by children is only possible under the supervision of adults.

Charging the battery

To charge the battery, insert the cable connector into the charging socket and connect the power supply unit to the mains. Charging can be done via a USB connection. To do this, disconnect the charger cable from the power supply and insert it into the USB connector of the included electronic device. In case of critical discharge and during charging, the mode indicators on the housing display flash. When the indicator stops flashing, it indicates that the battery has been charged. When using the device for the first time, the battery must be fully charged. A full charge is also required if the device has not been used for a long time. The battery takes 4 hours to fully charge.

Attention! You cannot use the device while the battery is charging.

Turning on the device

Open the tank cover, tilt the body so that the inlet is horizontal and perpendicular to the direction of the jet. Fill the tank with boiled, purified water. Close the tank cover. To turn on the device, press the power ON/OFF button once. To turn off the irrigator, press the ON/OFF button again.

Selecting the operation mode

The device can operate in three modes of different liquid pressure levels:

Normal (Normal)

Soft (Soft)

Massage (Massage)

The mode is selected by pressing the Mode button, which is confirmed by the inclusion of the corresponding light indicator on the handle of the device.

After switching on, the irrigator will operate in the same mode as it was last used before switching off the device.

Attention! At the beginning of using the device, it is recommended to use the soft mode. This will help you adapt to irrigation cleaning and avoid the risk of gum inflammation and bleeding.

Installation and replacement the nozzle

Install the nozzle in the connector located on the top of the irrigator and push it in until it clicks. To remove the nozzle, hold down the lock button and remove the nozzle from the device body.

Attention! Do not press the eject button while the device is in operation.

Irrigation mode

Position the nozzle at a 45-degree angle to the surface of the tooth. Turn on the irrigator, direct the water jet to the cleaning area. Water should flow naturally from the mouth to the sink. For best results, first treat the area of the side teeth, and then the front teeth. Move the nozzle along the contour of the gum, holding it in the interdental spaces.

Regardless of the mode, after operating for 120 seconds, the device automatically turns off.

After irrigation

Press the power off button. Drain the remaining liquid in the tank. If necessary, disconnect the tank from the appliance. To remove the tank, take the device body in one hand, and turn the tank counterclockwise with the other so that the marks on the lower part of the case and the upper part of the tank coincide.

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Attention! Some mouthwashes contain alcohol, which can damage the parts of the device for a long time and ultimately reduce its service life.

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Cleaning and maintenance

Do not immerse the device in water. It is forbidden to store liquid in the tank. Use a soft cloth and a soft, non-abrasive cleaner to clean the surface of the appliance. To facilitate the cleaning of the tank, it must be disconnected from the housing. With regular use of hard water, calcareous deposits may form on the parts of the device. To eliminate them, use a special Revyline product for cleaning the irrigator, according to the instructions for its use. The manufacturer recommends cleaning the irrigator from one to three times every three months, depending on the hardness of the water.

To avoid injuries or accidents related to the operation of the device, use only standard parts and accessories provided by the manufacturer.

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Model	Revyline RL 200	Revyline RL 200 XL
Parameter*		
Liquid pressure level	From 210 to 700 kPa	from 210 to 700 kPa
Tank volume	200 ml	300 ml
Working modes	3	3
Ripple frequency	1400 – 1800 per minute	1400 – 1800 per minute
Housing moisture protection	IPX7	IPX7
Supply	Lithium Battery 1400 mAh	Lithium Battery 1400 mAh
Time of full charge	4 hours	4 hours
Battery life	Till 500 cycles (charge/discharge)	Till 500 cycles (charge/discharge)
Battery life on a single charge	Minimum 60 minutes	Minimum 60 minutes
Operating time at full tank	At normal mode (Normal) 60 sec At soft mode (Soft) 80 sec At massage mode (Massage) 75 sec	At normal mode (Normal) 80 sec At soft mode (Soft) 120 sec At massage mode (Massage) 110 sec
Auto-Off	120 sec	120 sec
Power consumption, no more than	5 Vt	5 Vt
Noise level	≤ 78 dB	≤ 78 dB

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Bleeding during irrigation.

Solution: Bleeding is a characteristic symptom of gum inflammation. First of all, it is necessary to reduce the pressure level of the jet. The irrigator removes plaque and bacteria from the gums. If bleeding continues during further use, it is necessary to seek the advice of a dentist.

How to prevent water splashing?

Solution: During irrigation, do not stand upright with your mouth wide open. Lean towards the sink, open your mouth slightly. Water will naturally flow out of the mouth.

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	"wire-to-ground" scheme ± 2 kV		
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	The quality of electrical energy in the network is in accordance with the typical conditions of a commercial or hospital environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.
NOTE: U_H – the voltage level of the electrical network until the moment of applying the test action.			

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Note 1: For frequencies between 80 MHz and 800 MHz, a higher field strength value is applied.
Note 2: Expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people.

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证明书的网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

INSTRUCTION FOR USE

Irrigator oral dental portable model Revyline RL 800

Name:

Irrigator oral dental portable model Revyline RL 800 in composition:

1. main unit – 1 pc
2. water tank – 1 pc
3. standard nozzle – 3 pc
4. periodontal nozzle - 1 pc
5. orthodontic nozzle-1 pc
6. Tongue cleaning nozzle-1 pc
7. implants nozzle-1 pc
8. Handle – 1 pc.
9. instruction for use – 1 pc

Accessories:

1. case – 1 pc
2. USB adapter – 1 pc
3. USB cable – 1 pc

(hereinafter-irrigator, product, device, irrigator , Revyline RL 800)

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IMPORTANT!!!

Before using the Revyline RL 800 irrigator for the first time, read the operating instructions carefully and save them.

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1. Do not immerse the appliance in water or other liquid.
2. Do not use the irrigator while showering.
3. Avoid getting water in the charging port, getting wet, or damaging the USB cable.
4. Do not store the appliance in places where it may fall into a water container or sink.
5. Do not direct the water flow under the tongue, into the ears, nose, or other sensitive areas. The device generates a powerful pulsating water jet, which can disrupt the integrity of the skin and mucous membrane.
6. Use only branded attachments and accessories.
7. Do not insert foreign objects into the nozzle socket.
8. Charge the device away from hot surfaces.
9. Do not turn on the irrigator in a room where oxygen or aerosols have previously been used.
10. Do not fill the tank with rinses containing iodine, oil, bleaching ingredients, as well as rinses that allow natural precipitation and rinses that do not completely dissolve in water. They reduce the efficiency and service life of the device.
11. Fill the tank with water or a solution recommended by your dentist.
12. Before using the irrigator, remove all removable structures from the oral cavity.
13. Do not use the device if there are wounds in the mouth or on the tongue.

14. The use of the irrigator by children is possible from the age of 6 only under the supervision of adults.

Charging the battery

To charge the battery, insert the micro USB cable into the connector located on the back of the case. Connect the USB connector of the cable to a network adapter or a working electronic device equipped with a USB connector. While the battery is being charged, the indicator above the on/off button flashes. When the indicator stops flashing, it indicates that the battery has been charged. When using the device for the first time, the battery must be fully charged. The battery takes 4-5 hours to fully charge, depending on the power source.

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TURNING ON THE DEVICE

Install the handle nozzle in the socket, fix the handle in the recess of the fitting. Make sure that the pressure mode level is in the minimum position. Fill the tank with warm purified water to the required level, install the tank on the device body. Remove the handle, install the nozzle and point it into the oral cavity. To turn on the device, press the power button once. To turn off the irrigator, press this button again.

SELECTING THE OPERATING MODE

The liquid pressure level is set using the slider resistor key located on the left side of the housing. The rightmost position of the key corresponds to the minimum pressure level, the leftmost position corresponds to the maximum pressure level.

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Attention! At the beginning of the application of the device, it is recommended to use the device at a low pressure. This will help you adapt to irrigation cleaning and avoid the risk of gum inflammation and bleeding.

INSTALLATION AND REPLACEMENT THE NOZZLE

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AFTER IRRIGATION

Press the power off button. Drain the remaining liquid in the tank. If necessary, charge the device.

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	Date of manufacture (month, year)
	Caution! Special utilization required!
	Working part type BF
	Device class II
	Sign of product circulation on the territory of markets of participation of Custom union (in accordance with TR TS 020/2011)

Cleaning and maintenance

Do not immerse the device in water. It is forbidden to store liquid in the tank. Use a soft cloth and a soft, non-abrasive cleaner to clean the surface of the appliance. To facilitate the cleaning of the tank, it must be disconnected from the housing. With regular use of hard water, calcareous deposits may form on the parts of the device. To eliminate them, use a special Revyline cleanings agent to clean the irrigator, according to the instructions for use. The manufacturer recommends cleaning the irrigator from one to three times every three months, depending on the hardness of the water. To avoid injuries or accidents related to the operation of the device, use only standard parts and accessories provided by the manufacturer.

Transportation and storage of the irrigator

The Revyline RL 800 irrigator is convenient to use in road conditions, on trips. For the convenience of storage and transportation, the design of the device allows you to significantly reduce its dimensions. To do this, disconnect all parts of the device from the housing. Place the handle with the nozzle, the nozzles, and the USB cable in the case. Turn the device body over and insert it into the tank, and place the accessory case on top. Put the design in a travel bag, fasten the zipper. The irrigator is ready for transportation.

Attention! Before preparing the irrigator for the "road" form, it is necessary to remove the remaining water by wiping all the parts of the device with a dry cloth.

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Model	Revyline RL 800
Parameter*	
Liquid pressure level	От 140 до 800 кПа
Tank volume	500 ml
Working modes	5
Ripple frequency	till 1800 per minute
Housing moisture protection	IPX7
Supply	Lithium Battery 1400 mAh 1400 mAh
Time of full charge	4-5 hours (depending on the power supply)
Battery life	Till 500 cycles (charge/discharge)
Battery life on a single charge	Minimum 60 minutes
Operating time at full tank	1 mode 150 sec 2 mode 140 sec 3 mode 130 sec 4 mode 90 sec 5 mode 80 sec
Auto-Off	150 sec
Power consumption, no more than	5 Vt
Noise level	≤ 78 дБ

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Electromagnetic emission test	Compliance	Electromagnetic environment-instructions
Radio interference by CISPR 11	Group 1	Dental oral irrigators with accessories use radio frequency energy only to perform internal functions. The level of radio frequency interference emission is low and is not likely to cause disruption to nearby electronic equipment.
Radio interference by CISPR 11	Class B	
Harmonic components of the current according to IEC 61000-3-2	Class A	
Voltage fluctuations IEC 61000-3-3	Comply	

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Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Electrostatic discharges (ESR) according to IEC 61000-4-2	Contact discharge ± 6 kV	Contact discharge ± 6 kV	The floors should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, the relative humidity should be at least 30%.
	Air discharge ± 8 kV	Air discharge ± 8 kV	
Nanosecond pulse interference according to IEC 61000-4-4	For ± 2 kV power supply lines For input/ output lines ± 1 kV	For power supply lines of ± 2 kV	The quality of electrical energy in the network in accordance with the typical conditions of a commercial or hospital environment.
Microsecond high-energy pulse interference according to IEC 61000-4-5	When applying interference according to the "wire-to-wire" scheme ± 1 kV When applying interference	When applying interference according to the "wire-to-wire" scheme ± 1 kV	The quality of the power supply should be consistent with a standard commercial or hospital environment.

	according to the "wire-to-ground" scheme ± 2 kV		
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	$<5\% U_H$ (voltage drop $>95\% U_H$) for 0.5 period. $40\% U_H$ (voltage drop $>60\% U_H$) for 5 periods. $70\% U_H$ (voltage drop $>30\% U_H$) for 25 periods. $<5\% U_H$ ($>95\%$ voltage drop U_H) for 5 seconds	The quality of electrical energy in the network is in accordance with the typical conditions of a commercial or hospital environment. If the user of the products needs to ensure continuous operation in conditions of possible interruptions of the mains voltage, it is recommended that the product is powered from an uninterruptible power supply.
Industrial frequency magnetic field (50/60 Hz) according to IEC 61000-4-8	3 A/m 50 Hz/60Hz	3 A/m 50 Hz/60Hz	Industrial frequency magnetic field levels should be provided in accordance with typical commercial or hospital environments.
NOTE: U_H – the voltage level of the electrical network until the moment of applying the test action.			

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Interference resistance test	Test level according to IEC 60601	Compliance level	Electromagnetic environment-instructions
Conductive interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 Radio frequency electromagnetic field according to IEC 61000-4-3	3 V (mean-square value) in the band from 150 kHz to 80 MHz	3 V	The distance between the mobile radiotelephone communication systems in use and any element of the Irrigator, including cables, should not be less than the recommended spatial separation, which is calculated in accordance with the following expressions for the frequency of the transmitter. Recommended spatial spacing: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ (80 MHz to 800 MHz) $d=2.3\sqrt{P}$ (800 MHz to 2.5 GHz), where d is the recommended spatial spacing, m; P – the rated maximum output power of the transmitter, the maximum output power of the transmitter, W, set by the manufacturer. The field strength in the distribution of radio waves from stationary radio transmitters, according to the results of observations of the electromagnetic environment a) should be lower than the level of compliance in each of the frequencies b). Interference

			may occur in the vicinity of equipment marked with the sign: 
Note 1: For frequencies between 80 MHz and 800 MHz, a higher field strength value is applied.			
Note 2: Expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects, and people.			
a: The field strength during the propagation of radio waves from stationary radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and terrestrial stationary radio stations, amateur radio stations, AM and FM radio transmitters, television transmitters cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the Irrigator exceed the applicable compliance levels, the operation of the Irrigator should be monitored to verify their normal functioning. If the monitoring process reveals a deviation from normal operation, then additional measures may need to be taken, such as reorienting or moving the Irrigator. b: Outside the 150 kHz to 80 MHz band, the field strength must be less than 3 V / m			

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