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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01431757

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of WUXI HISKY MEDICAL TECHNOLOGIES CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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Authorized
Signature: Shen Xiaohua

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验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

 HISKY MED	Wuxi Hisky Medical Technologies Co., Ltd.	Document No. HME2018061
	Title of the document: User's Manual for a medical device "Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, model Pro 3000 advanced search, Pro 1 000"	Document vers 01
Date of change:	Date of entry into force: "23" November 2020	Page 1 of 13

User's Manual for a medical device "Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, model Pro 3000 advanced search, Pro 1000"

Approved

Director General

Wuxi Hisky Medical Technologies Co., Ltd.

Shao Jinhua

Shao Jinhua 2023.07.13
(signature)

Manufacturer:

Wuxi Hisky Medical Technologies Co., Ltd., China

Contents

Warning.....	VI
Chapter 1. Name and configuration of the medical device.....	1
Chapter 2. Information about the manufacturer of the medical device and its authorized representative	2
2.1 Manufacturer.....	2
2.2. Authorized representative of the manufacturer	2
Chapter 3. Purpose and potential users of the medical device.....	2
Chapter 4. Precautions for use	3
4.1 Security	3
4.1.1 Biosecurity	3
4.1.2 Indications and contraindications.....	3
4.1.2.1 Indications.....	3
4.1.2.2 Contraindications	3
4.1.3 Ultrasonic Safety.....	3
4.1.4 Data on medicinal products in the composition of a medical device.....	4
4.1.5 Electrical safety classification.....	4
4.1.6 EMC Declaration of Conformity	4
4.2 Safety regulations	4
4.3 Working conditions	7
4.4 Conditions of transportation and storage	8
4.5 Power supply.....	8
4.6 Precautions for use.....	9
4.7 Information about safety symbols.....	10
4.8 Terms of privacy and security	12
Chapter 5. Introduction to the system	13
5.1 Review	13
5.2 Intended use and indications.....	14
5.2.1 Intended use	14
5.2.2 Indications.....	14
5.3 Performance and features.....	14
5.3.1 Functional characteristics.....	14
5.3.2 Probe type	16
5.3.3 Display of elasticity test parameters	16
5.3.4 Probe pressure reading	16
5.3.5 Patient Information Management.....	16
5.3.6 Request function	17
5.3.7 Multi-user management	17
5.3.8 Automatic system diagnostic function	17
5.3.9 Regular compressed storage.....	17
5.3.10 Advanced compressed data storage.....	17
5.3.11 Smart cleaning function	17
5.3.12 HTTP and JSON support feature	17

5.3.13 Advanced search function	17
5.3.14 Visualization functions.....	17
5.3.15 Print function	18
5.3.16 Reporting function	18
5.3.17 Monitor	18
5.3.18 Main unit dimensions.....	18
5.3.19 Total weight of the product	19
5.3.21 General view of the product.....	19
5.4 Composition and components of the product.....	20
5.4.1 Main unit and component parts of the product	20
5.4.1.1. Main unit of the product.....	21
5.4.1.2 UAP Software	21
5.4.1.3 Imaging probe	22
5.4.1.4 Fibrosis Scanning Probe.....	23
5.4.1.5 Power cable 1.8 m.....	23
5.4.1.6 Foot switch.....	24
5.4.1.7 Electric fuse	24
5.4.1.8 Protective earth wire 2.2 m.	25
5.4.2 System components	26
5.4.3 Peripheral interface	27
5.4.4 Indicator light.....	27
5.5 Control panel.....	28
5.5.1 Function keys.....	28
5.5.2 Trackball	30
5.5.3 Monitor setup	30
5.5.3.1 Adjustment button	30
5.5.3.2 Adjusting of contrast and brightness.....	30
5.5.4 Power switch.....	31
5.5.5 Working part	31
5.6 Installation	31
5.6.1 Workspace.....	31
5.6.2 Equipotentiality Terminal.....	31
5.6.3 Probe connection.....	31
5.6.3.1 Connecting and disconnecting the probe for fibrosis scanning.....	32
5.6.3.2 Connecting and disconnecting an imaging probe	33
5.6.4 Foot switch connection	33
5.6.5 Connecting accessories	34
5.7 Movement and processing	34
5.7.1 Procedure	34
5.7.2 Additional Instructions.....	34
Chapter 6. System operation.....	36
6.1 Statement	36
6.2 Launch	36
6.2.1 Checks before use	36

6.2.2 Power on	36
6.3 Login.....	37
6.3.1 Login interface.....	37
6.3.2 Administrator login	37
6.3.3 Physician login.....	38
6.4 Interface operations.....	39
6.4.1 Create a new patient case	39
6.4.2 Case search.....	40
6.4.3 Deleting a case	42
6.4.4 Exporting Patient and Exam Data	42
6.4.5 Data interface	44
6.4.6 Data pinning.....	44
6.4.6.1 Configuration method	44
6.4.6.2 Fastening method	44
Chapter 7. Liver Elastography Detection.....	46
7.1 Preparing for detection.....	46
7.2 Ultrasound Imaging	46
7.2.1 Entering Imaging Mode	46
7.2.2 Imaging	46
7.2.3 Location defining for detection.....	48
7.2.3.1 Description of ultrasound images of the liver	48
7.2.3.2 Requirements for suitable localization.....	48
7.2.3.3 Image analysis of common errors	49
7.3 Elastography detection.....	51
7.3.1 Entering elastography mode.....	51
7.3.2 Detection launch	51
7.3.3 Mode M and mode A	53
7.3.4 Elasticity modulus.....	55
7.3.5 Detection results.....	57
7.3.5.1 Stiffness.....	57
7.3.5.2 Performance indicator	57
7.3.5.3 Ultrasonic attenuation parameter	58
7.3.5.4 Completion of detection.....	58
7.4 Printing a report	59
Chapter 8. Ultrasound Scan	59
8.1 Display modes.....	59
8.1.1 Mode B	59
8.1.2 Mode B/B.....	62
8.1.3 Mode 4B	63
8.1.4 Mode B+M.....	63
8.1.5 Mode M.....	64
8.1.6 Mode B+C.....	64
8.1.7 Mode B+D	65
8.1.8 Mode B+C+D	66

8.2 General measurements and calculations	66
8.2.1 B, 2B General measurements.....	66
8.2.1.1 Angle measurement.....	67
8.2.1.2 Volume measurement.....	68
8.2.2 General measurement V+M, M	69
8.2.2.1 Heart rate measurement	69
8.2.2.2 Valve speed measurement	70
8.2.3 General measurement PW.....	71
8.2.3.1 Time measurement	71
8.2.3.2 Heart rate measurement	72
8.2.3.3 Speed measurement.....	72
8.2.3.4 Acceleration measurement	72
8.2.3.5 General measurements	73
8.2.3.6 Doppler trace.....	74
8.3 Special measurements and calculations	74
8.3.1 Obstetric measurements	74
8.3.1.1 Measurement menu	74
8.3.1.2 Early OB Menu	75
8.3.1.3 OB measurement menu at later dates	76
8.3.2 Introduction to Obstetrics.....	77
8.3.2.1 Selecting a default database	78
8.3.3 Cardiac measurements	82
8.4 Ultrasound image parameter settings	82
8.5 Printing a report	83
8.5.1 Saving/accessing an image.....	83
8.5.2 Printing an image scan report.....	84
8.5.3 Saving a report	85
8.6 Working with shared buttons	86
Chapter 9. Work after use	92
9.1 Turning off the main unit	92
9.2 Power off.....	92
9.3 Checking of the technical condition.....	92
Chapter 10. Acoustic power.....	93
10.1 Problems related to biological effects	93
10.2 Reasonable use statement	93
10.3 ALARA principle (achieving the result with the lowest possible impact)	93
10.4 Explanation of MI/TI	94
10.4.1 Basic information about MI and TI.....	94
10.4.2 MI/TI display	95
10.5 Acoustic power management	96
10.6 Acoustic power limits	96
Chapter 11. System Maintenance.....	97
11.1 Probe Maintenance.....	97
11.1.1 Precautionary measures.....	97

11.1.2 Storage	97
11.1.3 Cleaning and disinfection.....	98
11.1.3.1 Cleaning and disinfecting the imaging Probe.....	99
11.1.3.2 Cleaning the Fibrosis Scan Probe.....	99
11.1.3.3 Cleaning the trackball.....	99
11.1.3.4 Cleaning of the control panel	99
11.1.3.5 Cleaning of the footswitch	99
11.1.3.6 Cable cleaning.....	100
11.1.3.7 Cleaning of the main unit.....	100
11.2 Maintenance of the main unit.....	100
11.2.1 Precautionary measures.....	100
11.2.2 Maintenance	100
11.2.3 Daily checklist.....	101
11.2.4 Weekly checklist	101
11.3 Security check	101
11.3.1 Security check	101
11.3.2 Troubleshooting	102
11.4 System repair.....	104
11.5 Disposal and environmental protection.....	104
11.6 Packing.....	105
11.7 Data on the labeling packaging of a medical device	105
11.7.1. Main unit marking.....	
11.7.1.1 Manufacturer's mark, model: Pro1000.....	
11.7.1.2. Manufacturer's mark, model: Pro3000 advanced research	
11.7.2. Marking for delivery to the Russian Federation.....	
11.7.3. Transport marking for delivery to the Russian Federation	
11.8. About symbols	
11.9 Warranty and service life	
Annex A: Acoustic power data	
Appendix B: EMC Declaration of Conformity	
Application C: List of international and national standards.....	

Warning

Before using the medical device, be sure to read the entire contents of this manual, including sections regarding machine installation, connection of peripherals, safety devices, functional applications, and maintenance methods.

Although the design takes into account the safety of the operator and the examined patient, for the reliability and safety of the system, attention must be paid to the following:

1. This medical product "Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, in versions Pro1000, Pro3000 advanced research" (hereinafter referred to as the system, device, product) should not be used by persons who are not fully qualified and certified medical professionals.
2. The system is classified as a Type BF, Class I application in accordance with the relevant safety standards. Please use the system properly in accordance with the information described in the Safety section, subsection 4.1 of this manual.
3. Do not make changes or modifications to the software or hardware of this system. If changes need to be made, please notify HISKY.
4. It is the customer's responsibility to maintain and manage the system after delivery.
5. The system has been debugged before delivery. Do not set up any preconfigured adjustable components.
6. In case of any abnormality during use, please turn off the power and inform the designated contact person.
7. The power cable of the system must be plugged into a grounded outlet. Do not remove or dismantle the protective earth wire.
8. If the system is to be electrically or mechanically connected to equipment from other suppliers, please check that such equipment complies with the current IEC60601-1 standard and re-check the leakage current and other safety features of the entire system to avoid any potential hazard resulting from the imposition of current leaks.
9. If the system is operated with any high-frequency surgical equipment, the operator must pay special attention, since the system is not equipped with special protective equipment.
10. You need to back up your important data to an external storage medium. Please back up important data on the system periodically, for example using a flash drive.
11. HISKY is not responsible for the loss of data stored in the memory of this system

caused by operator error or accidents.

12. A lithium battery is installed on the main board of the computer for the functioning of the clock memory. Please choose the same battery model for replacement. Do not immerse the lithium battery in water or place it at 100°C or higher at any time. Do not intentionally replace the lithium battery as it may be explosive. Dispose of the used battery in accordance with local rules or regulations.

13. This manual contains warnings about anticipated potential hazards, but you should always be aware of hazards other than those listed. HISKY is not liable for damage or loss resulting from negligence or disregard for the precautions and operating instructions described in this user manual.

14. If you change the manager of this system, be sure to hand over this manual to the new manager. If the user manual is not expected to be used for a long time, it should be retained for future reference.

15. The purpose of this system is to provide physicians with data for clinical diagnosis. The physician is responsible for diagnostic procedures. HISKY is not responsible for the results of diagnostic procedures.

 **Attention!** *We are not responsible for any damage and/or injury resulting from user actions that do not comply with the requirements of this manual or the instructions of the manufacturer/service center.*

User's Manual for a medical device " Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, in versions Pro1000, Pro3000 advanced research" System software and hardware are constantly updated, so the information contained here may differ slightly from your system.

Document No. HME2018061

Operating manual, version no. 01.

Release date November 23, 2020.

Chapter 1. Name and configuration of the medical device

" Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, in versions Pro1000, Pro3000 advanced research" (Further system, device, product).

Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis version Pro1000, consisting of:

1. Main unit - 1 pc.;
2. Monitor LCD 21.5" – 1 pc.;
3. UAP software -1pc;
4. Visualization probe YD-CX – 1 piece;
5. Fibrosis scanning probe XW-01 – 1 pc.;
6. Power cable 1.8 m - 1 pc.;
7. Foot switch, model TFS-201KD-2BWP-L - 1 pc.;
8. Fuse, brand T6.3AH250V - 2 pcs.;
9. Protective earth wire 2.2 m - 1 pc.;
10. Operation manual - 1 copy.

Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis version Pro3000 advanced research, consisting of:

1. Main unit - 1 pc.;
2. Monitor LCD 21.5" – 1 pc.;
3. UAP software -1pc;
4. Visualization probe YD-CX – 1 piece;
5. Fibrosis scanning probe XW-01 – 1 pc.;
6. Power cable 1.8 m - 1 pc.;
7. Foot switch, model TFS-201KD-2BWP-L - 1 pc.;
8. Fuse, brand T6.3AH250V - 2 pcs.;
9. Protective earth wire 2.2 m - 1 pc.;
10. Operation manual - 1 copy.

Description and (or) list of possible versions (configurations) of the considered medical product.

This medical device "Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, in versions Pro1000, Pro3000 advanced research" is available in 2 versions: model Pro1000 and model Pro3000 advanced research.

The Pro1000 model differs from the Pro3000 model in that with the Pro3000 model, the user can search the patient's medical record by entry period using the advanced search function, while the Pro1000 model does not have the advanced search function.

Chapter 2. Information about the manufacturer of the medical device and its authorized representative

2.1 Manufacturer

Name: Wuxi Hisky Medical Technologies Co., Ltd.

Legal address: Room B401, 530 Plaza, University Science Park, Taihu International Science & Technology Park, 214135 Wuxi, Jiangsu Province, PEOPLE'S REPUBLIC OF CHINA.

12.2. Authorized representative of the manufacturer

Limited Liability Company TESTSERVICE YUZHNY (LLC TESTSERVICE YUZHNY), Russia

Legal address: 143982, Moscow region, Balashikha, Yuzhnaya str. (Kuchino Mkr.), 9, floor 1, room. No 41. Telephone: +7 (499) 673-05-06. E-mail: regzn_0506@mail.ru

Chapter 3. Purpose and potential users of the medical device

The product is designed to measure liver stiffness/attenuation of ultrasonic waves in tissues based on transient elastography for the assessment/diagnosis of liver disease fibrosis/steatosis.

The Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, is designed to assess the density of liver tissue. Shear wave technology is used to diagnose liver fibrosis. The shear wave velocity is directly related to the elastic modulus of the tissue. Since fibrosis leads to hardening of tissues and, consequently, to a change in the elastic modulus, quantitative determination of the shear wave velocity allows an assessment of the level of fibrosis. The product also provides the imaging mode in order to guide appropriate placement of the shear wave assessment transducer.

The device is intended for use in medical and preventive medical institutions.

Potential users: Qualified and certified medical professionals with professional responsibilities that include performing ultrasound examinations in the field of gastroenterology.

Chapter 4. Precautions for use

4.1 Security

4.1.1 Biosecurity

The system is in contact with the human skin through an ultrasonic gel. For examination, the patient needs to apply an ultrasound gel. Therefore, it is necessary to use a qualified ultrasonic gel. Apply the ultrasound gel until the required clinical information is obtained.

4.1.2 Indications and contraindications

4.1.2.1 Indications

The device is designed for non-invasive measurement of the shear wave velocity of 50 Hz in the liver and assessment of the density, as well as the ultrasound attenuation parameter (UAP).

Shear wave velocity and density, as well as UAP, can be used as an aid in diagnosing and monitoring patients with liver disease as part of an overall liver assessment.

The device is designed for general-purpose pulse-echo ultrasound imaging and Doppler analysis of the human body. It can be used for the following areas of the body: Abdominal cavity, including the localization of the liver.

4.1.2.2 Contraindications

The device is not applicable to the wound at the measurement point.

Note:

1. Neonatal and pregnant patients should follow their physician's recommendations for a fibrosis screening test.
2. The presence of ascites between the transducer and the liver may interfere with the measurement.

4.1.3 Ultrasonic Safety

It is necessary to complete the detection as quickly as possible and use the minimum power prudently. If the scan is not completed, pause the probes or remove them from the patient. Do not keep the probes in the same position on the patient. The acoustic power level of the system is below the allowed value. The system is safe for users. See Appendix A for acoustic power ratings.

4.1.4 Data on medicinal products in the composition of a medical device

Medicinal products are not used as part of a medical device.

4.1.5 Electrical safety classification

1. The system meets the following security criteria:
 - 1) Working part type BF, class I;
 - 2) Continuous operation in operating mode (> 8 h);
 - 3) Start-up time of the product is less than one minute;
 - 4) Degree of protection of the main unit: IPX0 (enclosed equipment without liquid protection);
 - 5) Degree of protection of the probe for scanning fibrosis: IPX4.
 - 6) Degree of protection of the image probe: IPX4.
 - 7) Degree of protection of the foot switch IPX4;
 - 8) Not available in atmospheres containing a mixture of flammable anesthetic gas and air and/or oxygen and nitrous oxide.
2. Use of the system requires connection to a protective earth conductor. Make sure the system is properly grounded during use. Violation will be considered improper action and the user is responsible for its consequences, such as security risk and image violation.

4.1.6 EMC Declaration of Conformity

This system complies with the IEC60601-1-2 standard: 2014, and its RF emissions comply with CISPR11 A. See the Declaration of EMC Compliance in Appendix B.

4.2 Safety regulations

For the protection of patients and operators, please read the following safety instructions carefully before use.

⚠ Danger:

1. *DO NOT use flammable gases such as anesthetic gas or hydrogen or flammable liquids such as ethanol near this system as there is a risk of explosion.*
2. *The system is not intended for ophthalmic use or any of the possible ways in which the acoustic beam passes through the eye.*

⚠ Warning:

1. *The system is not equipped with a water protection element. Please do not place or touch the system where water can enter the system. Avoid any splash or water intrusion into the system. Otherwise, there is a risk of electric shock. In the event of accidental splashing, please contact our customer service department.*
2. *Correct power supply to ground according to the user manual. Otherwise, there is a risk of electric shock. It is forbidden to connect the protective earth wire to any gas or water pipe to prevent poor earthing or explosion hazard.*
3. *If any crack occurs on the surface of the product in contact with the human body, stop using such a probe immediately and contact our customer service department. Using a cracked probe will result in an electrical shock hazard.*
4. *Clean up any remaining chemical or gas after sanitizing accessories. Residual reagent or gas will not only damage accessories but may also cause injury.*
5. *Do not subject the system to excessive vibration (eg during transport). Otherwise, some mechanical parts such as the swivel wheel will be damaged. If the system needs to be moved regularly on uneven ground, please contact our customer service department.*
6. *Please do not use probes not supplied by us. If such a probe is used, the system and/or the probe may be damaged, even leading to fire and other accidents.*
7. *Do not open the case or panel. If the case is opened while the system is running, there is a danger of a short circuit or electric shock.*
8. *To avoid operator hand or wrist fatigue, please follow the instructions below: Avoid repeated operations for a long time and periodically rest or do other activities; do not use the system in an uncomfortable position. This will cause fatigue in the operator's wrists or hands.*
9.  This indicates that the system is classified as type BF equipment according to the classification of protection against electric shock. *Before using the system, check that all components are connected correctly and are working properly. Incorrect connection or abnormal operation may result in electric shock.*
10. *Do not allow the patient to touch live parts (such as signal input/output ports, etc.) as this may result in electric shock.*
11. *Don't hit the probe. Use of a defective probe may result in electric shock.*

12. *Since the use of defibrillators and high-frequency surgical equipment will cause interference with the instrument, it is not recommended to use the instrument in conjunction with such equipment.*
13. *Do not modify products, including system components, software, cables, etc. User-modified hardware can cause security issues and slow system performance.*
14. *The product must be disposed of in accordance with local disposal laws.*
15. *If the images are unstable or do not match, stop using the product. An intermittent scan sequence indicates hardware failures that must be corrected before use.*
16. *Only the special power cable supplied by HISKY may be used.*
17. *Warning: To avoid the risk of electric shock, this equipment must only be connected to a mains supply with a protective earth.*
18. *Connecting electrical equipment to the MSO creates a SYSTEM ME and the result can be a reduction in safety.*

**Attention:**

1. *Position the system horizontally and secure the swivel wheels to avoid impact injury to the patient or operator due to the system sliding.*
2. *The system can only move forward and backward. When moved to the left and right, the system may tip over.*
3. *Equipment connected to analog or digital interfaces must comply with the relevant IEC standards (e.g. IEC 60950 for data processing equipment and IEC 60601-1 for medical equipment). In addition, all configurations must comply with the current version of the IEC 60601-1 SYSTEMS standard. Anyone who connects additional equipment to the signal input part or the signal output part configures the medical system and is therefore responsible for ensuring that the system complies with the current version of the requirements of the IEC 60601-1 system standard. If in doubt, consult the technical service department or your local representative.*
4. *If the system is required to connect electronic or mechanical components manufactured by other companies, please consult us or our distributor in advance.*
5. *Don't place anything on the monitor. Otherwise, it will create an additional load on the system or lead to the danger of falling.*
6. *Always keep the system dry. Do not suddenly move the system from a relatively cold*

place to a warm place. Otherwise, it will cause condensation or water droplets to form, which may cause a short circuit.

7. If a circuit protector trips or a fuse blows, this indicates a system or peripheral failure. Please contact our sales office, customer service department or agent. Don't try to solve the problem yourself.

8. Conventional ultrasound examination will not lead to the risk of skin burns even if the transducer surface temperature is higher than the body temperature of the examined patient under various ambient temperatures and examination modes. However, do not place the probe on the same part of the patient's body for a long time to avoid skin burn. Minimize examination time as long as diagnostic conditions are met. Turn off the power to the system and contact us if there is a clear increase in the temperature of the probes.

9.  This indicates items that require attention. Please read the related information in the User Manual before using the system.

10. The working probes supplied with the system produce very low acoustic power and the mechanical and thermal effects are below the requirements of the applicable national standards ($MI < 1.0$; $TIS < 1.0$). Thus, they will not damage the fabric. However, please stop the probes after the examination to avoid unnecessary acoustic output.

11. The use of this system may cause the temperature in the narrow room to rise, please pay attention to the ventilation of the room.

12. If you need to dispose of the system or accessories, please contact HISKY or your local service representative. HISKY accepts no responsibility for any mishandling of this system without consultation that could result in loss or accident.

13. When using an ultrasonic probe, sterile gloves must be worn to prevent infection.

14. Please use the bonding agent in accordance with the regulations, and also use the ultrasonic bonding agent correctly.

15. Probes and other parts are not disinfected before leaving the factory, please treat them according to the instructions in the manual before use.

16. If you need to wear gloves when using the touch screen, please use gloves that can be used with capacitive touch.

4.3 Working conditions

1. Permissible normal operating conditions of the system:

- 1) Ambient temperature: $+5^{\circ}\text{C} \sim +40^{\circ}\text{C}$;
- 2) Relative humidity: 30%~80%;

- 3) Atmosphere pressure: 860hPa~1060hPa;
2. Strong sources of radiation or strong electromagnetic waves, such as electromagnetic waves from radio and television, may cause ghosting or noise. Thus, it is necessary to isolate the system from such sources of radiation or electromagnetic waves.
3. Do not place anything on the main unit.
4. When the system is not in use, turn off the power and cover the main unit with a dust cover.

1. *After transport, the temperature affected by the outdoor environment may be too high or too cold. Condensation may form due to temperature differences within the system. Before switching on, the system must adapt to the environment. If the outdoor temperature is lower than +10°C or higher than +40°C, the system must be left for some time in the working environment before use. If the temperature difference is 2.5°C, the storage time will be one hour. Repeat the procedure, if the temperature difference increases by 2.5°C, the equipment will not be able to work.*
2. *To avoid image interference, the working environment of the machine should be located away from the generator, X-ray machine, ultrasonic nebulizer, lithotripter, broadcasting station, television station, computer and transmission cables.*
3. *Ideally, the machine should be in an air-conditioned room.*

4.4 Conditions of transportation and storage

1. Transportation: The packaged system can be sent by regular transport. During transportation, protect the system from direct sunlight, rain, snow and mechanical damage.

Ambient temperature during transportation: -20°C +50°C.

Relative humidity during transportation: 30%-95%

2. Storage: The packaged system should be stored in a room with a temperature of -5°C to +40°C, relative humidity below 80%, and no corrosive substances.

4.5 Power supply

1. Power supply specification: 100-240V~, 50Hz/60Hz.
2. Main unit power: 350 VA
3. Main unit voltage: fluctuates within $\pm 10\%$ to prevent possible damage to the system.
4. Protective Grounding: Make sure the system is firmly connected to the power outlet during use.

⚠ Warning:

1. *Make sure the power supply meets the requirements. If you use an incorrect power supply, it may damage the product.*
2. *Depending on different regions, the system will have different power requirements. Please refer to the information on the back of the system for specific requirements.*

4.6 Precautions for use





1. Avoid using the system near a strong electromagnetic field such as a transformer. Otherwise, such conditions will adversely affect the monitor.
2. Avoid using the system near equipment with high frequency radiation, such as a cell phone. Otherwise, the system may fail or be adversely affected.
3. Avoid using the system in the following places to prevent damage to the system:
 - 1) Exposure to direct sunlight.
 - 2) Susceptibility to rapid temperature changes.
 - 3) Exposure to large amounts of dust.
 - 4) Frequent exposure to vibration.
 - 5) Placement near a heat source.
 - 6) Exposure to high humidity.
4. Do not restart the system immediately after shutting down, respect the mandatory minimum interval of 5 minutes.
5. After use, clean the ultrasonic gel on the probe and insert the probe into the probe holders. Otherwise, the moisture contained in the gel may enter the probes and adversely affect the performance and safety of the probes.
6. Insert and remove probes when the system is closed. Otherwise, it may cause problems.
7. The user can record registration information (including information about doctors, patients, etc.). To ensure data security, please back up your data on an external device. The data stored in the system may be lost in case of incorrect operation.
8. Transport precautions: Hold the handle when moving the system. If you hold on to other parts, the system will be damaged due to improper loading.
9. Avoid applying external force to the control panel, such as leaning on the panel. Otherwise, it may cause damage.
10. Take good care of the silicone membranes on your keyboard. Be careful to avoid scratches with sharp objects.

11. Peripheral power connectors on the main unit provide power to optional peripheral equipment. Do not connect other equipment to this connector. Otherwise, the system may fail if they exceed the output power of the power supply. The maximum output power is 100 VA for each of these connectors.
12. The maximum output power is 100 VA for each of these connectors. In this case, it is necessary to ensure good ventilation of the room.
13. A replacement fuse (T6.3AH250V) is installed in the system and must be replaced by professional service technicians or our authorized technicians.

4.7 Information about safety symbols

The system uses several of the symbols contained in the manual, including Danger, Warning, Attention, and Note. The meaning of these signs is described as follows:

Symbol	Value
 Hazard	Indicates imminent hazards that could result in death or serious injury.
 Warning	Indicates a potential hazard that could result in death or serious injury.
 Attention	Indicates a potential hazard that could result in minor or moderate injury.
Note	Indicates instructions that operators must follow to avoid damage to the system and to ensure normal and efficient operation and function of the system.
	Global warning
	Working part type BF
OFF.	OFF for NETWORK management only
ON.	ON for NETWORK management only
	Equipotentiality
	Attention: Electric Shock Hazard"
	Alternating current
IPX0	Not protected from water
IPX1	Protected against vertically falling water drops
IPX4	Splash-proof
IPX7	Protected against the effects of temporary immersion in water
IPX8	Protected against the effects of prolonged immersion in water

Symbol	Value
	Do not throw into trash. Use a separate collection for electrical and electronic equipment. Do not throw into trash. Use a separate
	See instruction manual/brochure NOTE on ME EQUIPMENT "Follow instructions for use"
	Manufacturer's declaration of conformity of the product with the applicable EEC directives, the European Notified Body is TÜV SÜD Product Service GmbH.
	Fuse
	Foot switch connector
	Top
	Fragile
	Protect from moisture
	Keep away from sunlight
	Stacking is prohibited;
	Ambient temperature for transportation: -20°C-50°C
	Ambient temperature for storage: -5°C-40°C
	Relative humidity during transportation: 30%-95%
	Relative humidity during storage: %
	Atmospheric pressure for transportation: 500hPa+1060hPa
	Atmospheric pressure for storage: 860hPa+1060hPa
	Main unit serial number
	Date of manufacture
	Manufacturer

Symbol	Value
	Authorized representative in European Community
	Additional information

4.8 Terms of privacy and security

1. Our service engineer can log into the background with a secret key to access the Windows desktop and operating system.
2. Operating systems are configured with security restrictions, including restricting Internet-connected programs from opening the firewall.
3. The system of non-essential, vulnerable network communication ports has been disabled.

Chapter 5. Introduction to the system

5.1 Review

The Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, is designed to assess liver tissue stiffness. Shear wave technology is used to diagnose liver fibrosis. The shear wave velocity is directly related to the elastic modulus of the tissue. Since fibrosis leads to hardening of tissues and, consequently, to a change in the elastic modulus, quantitative determination of the shear wave velocity allows an assessment of the level of fibrosis. The product also provides the imaging mode in order to guide appropriate placement of the shear wave assessment transducer.

The fibrosis scanning probe and the imaging probe are mainly composed of an ultrasonic lens, a housing and a cable. The acoustic lens is in contact with the surface of the human body through a bonding agent. The acoustic lens is made of silicone rubber, RTV630 brand, manufactured by Shenzhen Ultrasonic Electronics Co., Ltd, China. The surface of the fibrosis scanning probe is composed of synthetic resin (Epoxy polymer), grade 5010A/5010B+PF-017, manufactured by Wuxi Lanhui Ultrasonic Electronic Equipment Factory, China. Probe housings are made of ABS (Acrylonitrile Butadiene Styrene), grade PA-765A(+), manufactured by CHI MEI CORPORATION, China. Such materials are biologically compatible.

On fig. 5-1 there is a sketch of the main unit of the system.


Fig. 5-1 Main unit with 21.5" LCD monitor.

5.2 Intended use and Indications

5.2.1 Intended use

The Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis, a transient elastography product, is designed to measure shear wave velocity of 50 Hz and evaluate tissue stiffness and ultrasonic attenuation parameter (UAP) in internal body structures.

5.2.2 Indications

The Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis is designed to non-invasively measure the shear wave velocity of 50 Hz in the liver and evaluate the stiffness as well as the ultrasound attenuation parameter (UAP).

Shear wave velocity and density, as well as UAP, can be used as an aid in diagnosing and monitoring patients with liver disease as part of an overall liver assessment.

The Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis is designed for general purpose pulse echo ultrasound imaging and Doppler flow analysis of the human body. The system can be used in the following applications: Abdominal cavity, including localization of the liver.

5.3 Performance and features

5.3.1 Functional characteristics

1. Operating Parameters of Scanning Fibrosis Probe (XW-01)

Imaging Mode	Mode A Mode M Transition Elastography/Shear Wave
Source of ultrasound	Piezoelectric ultrasound source
Probe	Probe XW-01 (2.5 MHz) (single element ultrasonic transducer)
Elastography Mode	Transitional elastography
Source of mechanical vibration	External electromechanical vibrator
Determining Shear Wave Velocity	Post-processing
TE mode	Shear Wave Velocity and Tissue Density Measurements
TE display	Shear wave speed (0.8-5.2 m/s) Density (2.0-80 kPa) Interquartile range (IQR) and IQR/median ratio

Shear	(-2,3%) – (1,9%)
Accuracy	(0,0%) – (1,5%)
UAP display	UAP (90-450 dB/m) Interquartile range (IQR) and IQR/median ratio
Shear	(-1,5%) – (1,7%)
Accuracy	(0,1%) – (0,6%)

Note: Values obtained with phantoms 200GU-039-3, 201GU-039-3, 202GU-039-3, 3672.1-6-1, 940040-4665-1, 940100-4983-1, 940101-4986-1 and 940047- 4684-4 from CIRS et al.

2. Image/imaging probe operating parameters (convex, YD-CX)

No	Performance index	Requirements				
1	Probe frequency, MHz	3.5 (nominal)	2.0	2.5	3.0	5.0
2	Acoustic operating frequency, MHz	Deviation of the acoustic operating frequency from the nominal value within $\pm 15\%$.				
3	Examination depth, mm	≥ 160	≥ 180	≥ 180	≥ 180	≥ 100
4	Lateral resolution, mm	≤ 3 (depth ≤ 80) ≤ 4 ($80 < \text{depth} \leq 130$)	≤ 3 (depth ≤ 80) ≤ 4 ($80 < \text{depth} \leq 130$)	≤ 3 (depth ≤ 80) ≤ 4 ($80 < \text{depth} \leq 130$)	≤ 3 (depth ≤ 80) ≤ 4 ($80 < \text{depth} \leq 130$)	≤ 2 (depth ≤ 60)
5	Axial resolution, mm	≤ 2 (depth ≤ 80) ≤ 3 ($80 < \text{depth} \leq 130$)	≤ 2 (depth ≤ 130)	≤ 2 (depth ≤ 130)	≤ 2 (depth ≤ 130)	≤ 1 (depth ≤ 80)
6	Blind zone, mm	≤ 5	≤ 5	≤ 5	≤ 5	≤ 4
7	Cut thickness, mm	≤ 8	≤ 8	≤ 8	≤ 8	≤ 7
8	Accuracy of horizontal geometric position, %	≤ 15	≤ 15	≤ 15	≤ 15	≤ 15
9	Accuracy of longitudinal geometric position, %	≤ 10	≤ 10	≤ 10	≤ 10	≤ 10
10	Measurement error of circumference and area, %	± 20	± 20	± 20	± 20	± 20
11	Time display error in M mode, %	≤ 5	≤ 5	≤ 5	≤ 5	≤ 5

CF mode:

No	Performance	Unit	Requirements	
1	Probe frequency	MHz	2.0	3.0
2	Examination depth	mm	≥100	≥80

PW mode:

No	Performance	Unit	Requirements	
1	Probe frequency	MHz	2.0	3.0
2	Examination depth	mm	≥100	≥80

5.3.2 Probe type

1. Fibrosis Scanning Probe: Used in Elastography mode to examine liver elasticity to help scan liver fibrosis. It is applied to the surface of the patient's abdomen on the right.
2. Imaging probe: used in the abdominal imaging mode for universal pulse echo ultrasound imaging and Doppler flow analysis of the human body. Ultrasound imaging also helps to find a suitable site for temporal elastography examination.

5.3.3 Display of elasticity test parameters

The interface displays information about the examined patient, stiffness (kPa), quartile deviation, relative deviation (%), performance indicator (%), UAP (dB/m), number of examinations, number of successful examinations and pressure indicator.

5.3.4 Probe pressure reading

The transducer pressure indicator bar will turn green if the pressure is within normal limits during the test. Otherwise, it will turn red.

5.3.5 Patient Information Management

Create, modify, save, query, delete or otherwise process patient information.

5.3.6 Request function

You can search for a patient in the “Search”/“Patient” mode by keyword: ID, Name, Age, Gender, Phone or Time.

5.3.7 Multi-user management

You can set the resolution for multiple users.

5.3.8 Automatic system diagnostic function

The system can be detected and repaired automatically through the interface operation.

5.3.9 Regular compressed storage

The system can store more patient data after normal compression.

5.3.10 Advanced compressed data storage

The system can store more patient data after advanced compression.

5.3.11 Smart cleaning function

When performing Smart cleaning, if the storage capacity is below the threshold, the original data will be automatically deleted from the beginning.

5.3.12 HTTP and JSON support feature

The system can be connected to the server via HTTP protocol in JSON format.

5.3.13 Advanced search function

Searching for a patient by entered time period is only applicable to the Pro3000 model. The Pro1000 model differs from the Pro3000 model in that with the Pro3000 model, the user can search the patient's medical record by entry period using the advanced search function, while the Pro1000 model does not have the advanced search function.

5.3.14 Visualization functions

The imaging mode functions are shown in the table.

Imaging modes	A, B, B/B, 4B, M, B+M, CF, PW, Harmonic Image, Sync Display, Zoom, Depth, Beamline
---------------	--

Picture control features	Settings	Frequency, Focal Number, Dynamic Range, TGC, Gain, Acoustic Power, Harmonic Image, Center Line, Depth, Edge Enhancement, Up/Down, Left/Right	
	Setting parameters in CF mode	Pulse Repetition Frequency (PRF), Wall Filter, Packet Size, Color Level, CF Block Position, Acoustic Power, Baseline, Line Density, CF Block	
	Setting parameters in PW mode	Pulse Repetition Frequency (PRF), Display Sync, Wall Filter, Baseline, Sample Size, Acoustic Power, Scan frequency	
	Basic measurement functions	Perimeter	Distance
		Area	Length and width, ellipse
		Measuring in M mode	Heart rate, Time, Distance, Speed
		Volume, Angle, Diameter Reduction, Area Reduction, Tracking	
	Measurement in PW mode	Speed, Time, Acceleration, Heart Rate, General, Tracking	
	Software package	Abdomen, Obstetrics, Pelvis, Urology, Small organs, Vessels, Heart, Musculoskeletal system	
	Image Format	BMP, JPEG, AVI, DCM, PNG	
	Language	Chinese, English	
	Ports	USB, Ethernet, footswitch	
Administration		New Patient, Patient Archive, System Setup, Body Drawing, Print, Storage, Report, Patient Data Export, Search, Archive, etc.	

5.3.15 Print function

The print function includes a liver elasticity measurement report.

5.3.16 Reporting function

1. Access to image file, patient information and examination report from hard drive;
2. Access image file, patient information and exam report via USB stick.

5.3.17 Monitor

LCD display: 21.5" TFT LCD, weight 4,7 kg $\pm$ 10%.

5.3.18 Main unit dimensions

L $\times$ W $\times$ H=95 cm $\times$ 65 cm $\times$ 145 cm.

5.3.19 Total weight of the product

Total mass of the main unit: 81 kg $\pm$ 10%.

5.3.20 Effort of movement of the main block:

Effort of movement of the main block:

- with the brakes released, not more than 50 N.
- with the brakes on, no more than 295 N.

5.3.21 General view of the product


Fig. 5-2 General view of the apparatus.

5.4 Composition and components of the product

5.4.1 Main unit and component parts of the product

Name:	Quantity	Use
Main unit	1	Liver stiffness calculation using obtained data
Monitor LCD 21.5"	1	Displaying the control interface and scan data
UAP Software	1	Embedded software for ultrasound diagnostics
YD-CX Imaging Probe	1	Determination of the testing area
XW-01 Fibrosis Scanning Probe (2.5 MHz)	1	Liver stiffness test and UAP

Power cable 1.8 m.	1	Main unit power
Foot switch, model TFS-201KD-2BWP-L	1	Start of the survey
Fuse, brand T6.3AH250V	2	Overcurrent protection
Protective earth wire 2.2 m.	1	System grounding

5.4.1.1. Main unit of the product

Fig.5-3. The main unit of a medical device with a 21.5" LCD monitor.

5.4.1.2 UAP Software

The software for ultrasound diagnostics is a program designed to be installed in an application system based on such hardware as a database and a fibrosis scanning probe. Accordingly, for the implementation of the project, it is necessary to achieve the following goals at the system and management levels.

(1) Software system name: Pro3000/ Pro1000 software system.

The medical diagnostic system must meet the standards of state, provincial and municipal hospitals at various levels, and provide fast and accurate measurement and analysis of data, meet the requirements of auxiliary diagnostics, and provide effective assistance, so that the attending physician can make a timely diagnosis.

(2) The system should standardize and streamline the processes of reporting and summarizing patient information, maintain concentration and communication of patient information, improve management efficiency and quality, and improve the level of hospital management.

(3) The system should support the gradual implementation of the unification and standardization of rules and regulations.

(4) Software risk class:

The software of the device corresponds to the class - 2a of the potential risk of using a medical device.

(5) Software safety class:

In accordance with the hazard analysis of the medical device, the software system was assigned a safety class - A. Class A: No injury or damage to health is possible, which corresponds to the possible effects on the patient, user or other persons.

Key software definitions.

C#: a high-level programming language released by Microsoft, object-oriented and running on the NET Framework;

VS2010: Used by Microsoft Corp development environment, which is currently the most popular Windows platform application development environment;

CPWD: color Doppler analysis of the flow spectrum;

PWD: Pulsed Doppler;

CFM: Color Doppler Imaging.

5.4.1.3 Imaging probe

Used in the belly imaging mode for general purpose Pulse Echo Ultrasound and Doppler analysis of the human body. Ultrasound imaging also helps to find a suitable site for shear wave elastography.

Modes: A, B, B/B, 4B, M, B+M, CF, PW.


Fig.5-4. Imaging probe.

5.4.1.4 Fibrosis Scanning Probe

Used in elastography mode to examine liver elasticity to assist in scanning for liver fibrosis. It is applied to the surface of the patient's abdomen on the right.

Modes: A, M and elastography.


Fig.5-5. Fibrosis Scanning Probe

5.4.1.5 Power cable 1.8 m.

Designed to supply power to the main unit.


Fig.5-6. Power cable 1.8 m.

5.4.1.6 Foot switch

Designed to run a survey


Fig.5-7. Foot switch.

5.4.1.7 Electric fuse

Designed for surge protection.


Fig.5-8. Fuse, brand T6.3AH250V

5.4.1.8 Protective earth wire 2.2 m.

Designed to connect the device to the ground


Fig.5-9. Protective earth wire 2.2 m.

5.4.2 System components


Fig. 5-10 Front view


Fig. 5-11 Back view.

5.4.3 Peripheral interface

Peripheral interfaces see fig. 5-11, 5-12. Rear view including the following:

1. Power input port
2. Peripheral power port
3. USB ports
4. Ethernet port
5. Foot switch port
6. Fuse holder
7. Equipotentiality terminal
8. Main power switch


Fig.5-12. Photo of the back panel of the device.

5.4.4 Indicator light

1. Power button: Green light means the device is on.
2. LCD indicator: Green light means the LCD is working, red light means there is no signal on the LCD.
3. Keyboard indicator: The blue indicator means that the function key is enabled.

5.5 Control panel

The diagnostic system is designed to create a user-friendly clinical environment that is comfortable to work with. The system aims to minimize unnecessary operations of operators and provide convenient use.

5.5.1 Function keys

The location of the keyboard panel, panel color, key name, etc. may vary slightly depending on the model, but the corresponding function is the same.



Image	Name:	Function Introduction
	System Button [On/Off]	Press the key to turn on the system after turning on the main power; Press the key to turn off the system after use.
	PC keyboard	Used to enter document control characters.
	TGC adjustment slider	Provides 8-step local image gain control in real time. Gain increases as you move the slider from left to right.
	[Setup]	Switch from depth, focus and zoom to adjust the corresponding setting.
	[Menu]	Selecting a menu item.
	[B]	Enter imaging mode and gain adjustment.
	[Patient]	Creation of information about the patient.
	[Search]	Request for patient information.





Image	Name:	Function Introduction
	[Mode]	Switch between elastography scan and image scan.
	[Pause]	Enable or disable image scanning in imaging mode.
	[Snapshot]	Clipping video images in render mode.
	[Print]	Print an Elastography/Image Scan Report.
	[Cancel]	Cancel the current operation.
	[Save]	Image saving elastography scan and image scan.
	[probe]	Search the imaging probe in use.
	[Up]	View the last snapshot of an image in imaging mode.
	[Down]	View the next image snapshot in imaging mode.
	[Clear]	Clearing patient information.
	[Measurement]	Measurement of tissue dimensions after a pause in the imaging mode.
	[Input]	Left mouse button used with a trackball. If there is no special declaration, the [Enter] key mentioned in the manual refers to this key.
	[Exit]	The right mouse button used with a trackball.
	[F1] - [F10]	Function selection keys.

Image	Name:	Function* Introduction
	[B]	Enter B mode.
	[B/M]	Enter B M mode.
	[PW]	Enter P / mode.
	[CFM]	Enter C M mode.
	[CINE]	Entrance to the CINE overview (X-ray cinematography).

5.5.2 Trackball

Used according to the working status an menu to carry out all kinds of operating functions.

Management of a multiplex electronic ruler, position indicator, performance of measuring functions.

5.5.3 Monitor setup

5.5.3.1 Adjustment button

1. "Power" - turn on the monitor.
2. "Down" – down button.
3. "Up" - moving up button.
4. "Menu" – menu button.
5. "Auto" – automatic adjustment button

5.5.3.2 Adjusting of contrast and brightness

1. Automatic adjustment

After pressing the "Auto" button, the system will set a relatively better state of contrast and brightness.

2. Manual adjustment

- 1) Press the "menu" button to display the LCD monitor setting menu.
- 2) Press "Up" or "Down" to select brightness (or contrast) and press "Menu" to confirm your selection.
- 3) Press "Up" or "Down", to adjust the brightness (or contrast).
- 4) After reaching the desired brightness (or contrast), press the Menu button to confirm your selection.

- 5) Press “Up” or “Down” to select EIT menu, and press “Menu” to exit monitor settings.

5.5.4 Power switch

1. Main Power Switch: Located on the back of the device, it is used to turn the system power on or off.
2. System on/off button: on the control panel. Used to start or stop the system.

5.5.5 Working part

There are two working parts of the BF in the system: the imaging probe and the fibrosis scanning probe.

5.6 Installation

5.6.1 Workspace

Leave a space of at least 20 cm to the left and behind the system.

⚠ Attention! *Leave enough space to the left and behind the system. Otherwise, an increase in temperature in the system may lead to failure.*

5.6.2 Equipotentiality Terminal

The system power cable is a three-core cable. The equipotential terminal must be connected to the earth protection mechanism of the power supply. Please make sure the power supply ground protection is working properly.

⚠ Warning:

1. *Do not connect the system's three-wire power cord to any two-pole power outlet. Otherwise, there is a risk of electric shock*
2. *High input power, pay attention to operation safety.*

5.6.3 Probe connection

⚠ Attention!

1. *Connect and disconnect probes only in the off state. Otherwise, a failure may occur.*
2. *When connecting or disconnecting probes, place the probes in an appropriate holder and hang the probe cable in a cabinet rack to protect the probes from damage if dropped.*
3. *Do not use force to connect or disconnect the probe.*
4. *Always protect probes from impact, drops, scratches, damage, or other hazards.*

5. Use only approved contact fluid.
6. Do not remove the plug from the socket while the machine is in operation.
7. Do not pull or bend probe cables or other power cables.
8. Do not allow probes to come into contact with paint thinner, ethylene oxide, or other organic solvents.

Note: The use of our original probes is mandatory. Connecting other probes may damage the system and probes.

5.6.3.1 Connecting and disconnecting the probe for fibrosis scanning

The fibrosis scanning probe uses a self-locking plug-in connection. It can effectively resist vibration and shock, avoid disconnection when pulling the cable, ensure safe and reliable use, and facilitate the connection and disconnection of the probe in a very limited space.

1. Probe connection

⚠ Warning: In advance, check the probes cable and connectors for faults (e.g. cracks and disconnection). Use of defective probes will create a risk of electrical shock.

As shown in the following fig. 5-13, the right half is the plug that connects the fibrosis scan probe, and the left half is the socket permanently attached to the system case. When connecting the fibrosis probe, point the protrusion on the plug (indicated by the red dot) into the socket, and then insert the plug into the socket in the direction of the arrow and lock it.


Fig. 5-13 Probe connection

2. Disconnecting the probe

⚠ Attention! To unlock the probe, pull the release mechanism on the plug (as shown in Figure 5-14). Otherwise, pulling on a cable or other part will not disconnect the transducer and may damage the fibrosis probe.

If it is necessary to disconnect the probe, pull the release mechanism in the direction shown in fig. 5-7. Unlock the lock pin and remove the plug from the outlet.


Fig. 5-14 Disconnecting the probe

5.6.3.2 Connecting and disconnecting an imaging probe

1. Probe connection

⚠ Warning: In advance, check the probes cable and connectors for faults (e.g. cracks and disconnection). Use of defective probe will create a risk of electrical shock.

- 1) Release the probe connector lock switch. Hold the cable up and insert the plug vertically into the socket.
- 2) Rotate the probe connector lock switch 90° clockwise.
- 3) Check if the probe is blocked.

2. Disconnecting the probe

Rotate the probe connector lock switch 90° counterclockwise and pull the connector vertically.


Fig. 5-15 Probe connection Fig. 5-16. Probe blocking Fig. 5-17. Disabling the probe

5.6.4 Foot switch connection

The footswitch uses the same self-locking plug connection as the fibrosis probe. The red dot will help install and unlock the mechanism. Its connection and disconnection can be done according to the connection and disconnection method of the fibrosis scan probe.

5.6.5 Connecting accessories

1. Place the graphic printer in a stable position and connect the printer to the data port and graphic printer port on the back of the system with a printer cable.
2. Connect the graphics printer's power cable to a power source.
3. The monitor's articulated arm can be rotated 90 degrees clockwise or 90 degrees counterclockwise. The rotation range is not more than 180 degrees.

5.7 Movement and processing

Attention:

1. Before moving the system, check the following points:
 - 1) Turn off the power switch, unplug the power plug, and hang the power cord from the power cord hook on the back of the system.
 - 2) Place all probes in the holders and hang the probe cable in the rack.
 - 3) Disconnect and remove the footswitch.
 - 4) Set the monitor to its normal position.
 - 5) Start moving the system slowly and check if there are any deviations. In case of any deviations, please contact our sales office, customer service department or agent.
2. The system can only move forward and backward. If the system moves left and right, it may roll over.
3. The air circuit breaker installed on the back of the product is used to turn off the power. Please leave enough space for the product to be able to control the air circuit breaker.

5.7.1 Procedure

1. Unlock locked swivel wheels.
2. Hold the handle at the back of the system and push the system back and forth.
3. After moving, place the system in a stable position and lock the swivel wheels.

Note: After moving, check the system as it was before using it.


Fig. 5-18 Swivel wheel lock mechanism

5.7.2 Additional Instructions

If the system is moved to a remote location, follow these additional instructions:

1. Pack the system properly, preferably using the original packing materials, to avoid

shock, vibration, and extrusion during shipping.

2. Do not turn the system upside down during transport.

Chapter 6. System operation

6.1 Statement

Warning:

1. *The system must be operated by a trained medical professional.*
2. *Before use, read the system safety information described in Chapter One.*

6.2 Launch

 **Warning:** *Scheduled maintenance and inspection ensure safe and efficient use of the system. The following checks must be carried out before each start. If any abnormality is found, stop the system immediately and contact our sales office, customer service department or agent. Using a faulty system can harm the patient.*

6.2.1 Checks before use

Before turning on the power supply, perform the following checks:

1. Check all power supplies and connecting cables to make sure they are free of damage such as scratches and breaks.
2. Check the control panel, monitor, and base for any abnormalities such as cracks.
3. Check the probes and connections and make sure that the ultrasonic transmitting transducer does not have any abnormalities such as cracks and scratches, and that there is no loose or broken cable.

6.2.2 Power on

1. Please observe the following system launch order:
 - 1) Turn on the main power switch.
 - 2) Turn on the monitor.
 - 3) Press the system on/off button.
 2. Displays the system launch and initialization screen, followed by menu and image display. Check if the system starts normally.
 3. Check the surface of the working probe and make sure it is not overheated.
-  **Warning:** *Using an overheated probe may cause patient burns.*
4. Check the image to make sure there are no abnormalities such as abnormal noise or flicker.
 5. Check the control panel and make sure all buttons and knobs are working properly.

Warning: After using for one day, turn off the main power switch on the back of the product.

Warning: Deviation indicates defects in the system. In case of any deviations, please contact our sales office, customer service department or agent.

6.3 Login

6.3.1 Login interface

After starting, the login interface is shown in Fig. 6-1. Enter the relevant information in the username and password fields and click the "Login" button to log in.


Fig. 6-1 Login interface

For security reasons, users are divided into two groups: doctors and administrators. Different groups have different powers.

Additional information

1. The administrator assigns the username and password required to log in to the authorized physician. Enter the username and password to log in, which has a password of 8-24 characters, consisting of numbers and letters.
2. If the password is entered incorrectly more than 5 times during login, the user will be prevented from trying to log in again. The user will need to wait 10 minutes or reboot the device before logging in.

6.3.2 Administrator login

The administrator account is used to manage physician accounts and to add or remove any physician account. The login interface is shown in the following fig. 6-2.

Fig. 6-2 Administrator interface

After logging in, the administrator enters the management interface, as shown in Fig. 6-2. Here the administrator can create, modify or delete any user.

1. **User creation.** Click the "Create" button in the upper right corner. The control interface is changed to the interface shown in Fig. 6-3. In the "User Information" field, enter the user, password, name, gender, department, and other information. Click the "Save" button to save the entry.

Fig. 6-3 Creating a user

2. **User change.** Double-click the user you want to change in the user list. The information about the selected user is then displayed in the User Management field. Click the "Edit" button in the upper right corner. The control interface is changed to the interface shown in Fig. 6-4. After making changes, click the "Save" button to save the changes.

Fig.6-4 User change

3. **Deleting a user.** Double-click the user you want to delete in the user list. The information about the selected user is then displayed in the User Management field. Click the "Delete" button in the upper right corner. The control interface is changed to the interface shown in Fig. 6-5. After making changes, press the "Y (Enter)" button in the prompt dialog box to save the deletion.


Fig. 6-5 Deleting a user

6.3.3 Physician login

The doctor account is provided to the doctor who manages the system. After logging in through the doctor account, the doctor can see the interface as shown in the following fig. 6-6. Here, the doctor can perform fibrosis scans, ultrasound imaging, and other surgeries. After entering the interface, two options will be displayed: fibrosis scan and

image scan. Now press “Patient”, “Search” and other keys in the menu to control the system.


Fig. 6-6 Doctor's working interface

Additional information: *After logging in, the probe initialization process will occur, please wait for this process to complete, and then start the system. You will hear sounds during this process, this is normal.*

6.4 Interface operations

6.4.1 Create a new patient case

After logging in, the doctor enters the operation interface, as shown in Fig. 6-6. Here, click on the “Patient” symbol in the upper right corner of the screen or the function key [Patient]. Then the control interface goes to the patient case download interface, as shown in Fig. 6-7.

Fig. 6-7 Patient Button

After entering the patient information, press the “Enter” button directly and enter the image scanning interface. If the entry is incorrect, click the change icon “” and change the patient information. Then press the Enter button again.

Additional information:

1. Considering the circumstances in which the length of patient information characters may vary when entered, we have deliberately set the number of characters associated with new patient information to 100 bits in order to avoid problems for users due to the lack of character bits in actual use.
2. If there is a medical history in the system, you can double-click the existing medical history and download the information directly.
3. After successfully creating a new case, entering the default imaging interface, press the [Mode] button on the control panel to enter the fiber scanning interface.
4. When entering information about a case, you must fill in the menu spaces marked with an “*”, otherwise you will not be able to create a case. The "New Patient Case" menu is shown in fig. 6-8.


Fig. 6-8 New Patient Case Menu

Description of the "New Patient Case" menu:

"Clear": Clear information when entering a blank field.

"Enter": Confirm and save the new case information.

"Exit": exit the current state, return to the main screen.

6.4.2 Case search

The operation is used to search for a patient's medical history and medical results for comparison and analysis.

After logging in with a doctor account, enter the doctor working interface. If you want to search the history of patient information, please click "Search" in the upper right corner or hotkey "Search".

Fig. 6-9 Search button

1. Enter the relevant request information and click the Save button.
2. Information about medical records can be presented in two forms: an array of cards and an information list. Users can press the "Map" button to select the map array form, or "List" to select the information list form according to their preferences. The map array interface is shown in Fig. 6-10, the information list is shown in fig. 6-11.


Fig. 6-10 Main array interface


Fig. 6-11 Information list interface

- The medical record information supports the record appearing function. Click the "Merge by Patient" button to display the combined case data of the same patient for easy viewing. The "Combine by Patient" option is shown in fig. 6-12.

Fig. 6-12 "Merge by Patient" option

- The user can double-click the desired entries and view the diagnostic history.

6.4.3 Deleting a case

For the selected patient record, you can click "Delete" on the interface to delete the record for that time. The "Delete" button is shown in fig. 6-13.

Fig. 6-13 Delete Button

6.4.4 Exporting Patient and Exam Data

Additional information *The exported data does not contain sensitive information about patients or clients.*

- Click the "Export" button in the search interface to save the basic test data. The Export button is shown in fig. 6-14.

Fig. 6-14 Export button

- Choose the export path, you can choose local disk, you can also choose external U disk or external mobile hard disk. As shown in fig. 6-15.


Fig. 6-15 Export path

- Specify the name of the exported file as shown in Figure 6-16. The exported data is text in CSV format, which can be saved to a local drive, an external USB drive, or an external mobile hard drive.


Fig. 6-16 Export file name

- Click the "Save" button to save the exported data. The "Save" button is shown in fig. 6-17.


Fig. 6-17 Save Button

6.4.5 Data interface

1. Connect USB disk, mobile hard disk and other common external storage devices via USB interface for data exchange. Stored data is patient examination data and basic patient information.
Save formats: JPG, CSV. .
2. Connect to network via RJ45 network interface, upload survey report to SQL Server database, Oracle database and DB2 database server, download survey report via DICOM protocol, download survey report via HTTP protocol, and download survey report via SFTP protocol.

6.4.6 Data pinning

Users who need to upload and download patient data must set up the data pinning feature.

6.4.6.1 Configuration method

Please contact our after-sales service personnel to determine the data table information, such as the database format, and provide us with server information. We change the data pinning parameter in the configuration interface, at the same time in the pinning function configuration, according to the actual needs of uploading and downloading the function configuration.

6.4.6.2 Fastening method

1. In the new patient interface (as shown in Figure 6-7), enter the patient ID in the patient information input field, press the enter key on the PC keyboard, or check the box to the right of the down arrow icon  (see Figure 6-18), the system will automatically search for patient information on the hospital server and return to display in the current new patient interface corresponding to the input field.

Additional information: support database, DICOM worklist, HTTP and JSON.

***Required**

ID*	3
Name*	Amy
D.O.B.*	25 1995/6/19
Gender*	F
Height*	165 cm
Weight*	50 kg
Phone	
Clinic	
Physician	
Note	

Clear New Esc

Fig. 6-18 New patient interface

2. After the examination, the patient enters the report printing interface and clicks the download arrow, causing the diagnostic system to download the examination result or the report of the current patient to the hospital server;

Next (Enter)	Return (Esc)
	 

Fig. 6-19 Report Interface

Additional information: database, DICOM, SFTP, HTTP and JSON support.

Chapter 7. Liver Elastography Detection

7.1 Preparing for detection

Liver elastography detection by this instrument is applied to the right lobe of the liver. The patient lies on the bed during the operation of the scan procedure. The necessary positions of the patient are given below for better detection by the doctor:

1. Lying on the back with the right arm raised; expansion of the intercostal space if possible.
2. The patient continues to breathe evenly and does not speak for the time of detection.

For patients with narrow intercostal space:

1. Lying position, the body is bent to the left in type C.
2. The patient does not move the hip, the right leg is placed on the left.

The use of the ultrasound imaging function can provide an accurate selection of areas for detection, improving the reliability of elastography detection. You can also perform elastography detection directly without ultrasound imaging.

7.2 Ultrasound Imaging

7.2.1 Entering Imaging Mode

Perform an ultrasound scan before a liver fibrosis scan to more accurately position the liver tissue, select areas for examination, and get better results. After creating a new case, or press the hot key “Mode” to enter the imaging mode.

7.2.2 Imaging

Apply the gel to the transmitting surface of the ultrasound imaging probe and place the probe close to the skin over the liver.

Procedures:

1. Hold the probe with your right hand with your thumb on one side and the other four fingers on the other. The handle of the probe should be located closer to the area between the thumb and forefinger. The direction of the thumb should be the same as that of the inflection point.


Fig. 7-1 Probe hold

2. Touch the skin, starting from the xiphoid process, move down and to the right, you will find the 7th-9th intercostal spaces, as shown in Fig. 7-2.


Fig. 7-2, 7-9th intercostal spaces

3. Hold the probe and move along the intercostal spaces within the defined area. Usually the liver is located 2.5 cm below the midaxillary line and the 7th or 8th intercostal space, the thickness, specific location and shape of the liver will vary from person to person.

7.2.3 Location defining for detection

7.2.3.1 Description of ultrasound images of the liver

Fig. 7-3 Ultrasound images of the liver

The liver is an almost homogeneous tissue. An ultrasound image of the liver is shown in Fig. 7-3, where the fabric is marked with numbers:

1. Fibrous membrane of the liver
2. Vessels
3. Gallbladder

7.2.3.2 Requirements for suitable localization

Look at the outline of the fibrous tunic of the liver to make sure the liver is scanned. Within the white line near the center line of the image there is the fibrosis scan area. A suitable detection localization is shown in Fig. 7-4, and the specific requirements are as follows:

1. In the area of fibrosis scan there is a homogeneous liver parenchyma.
2. Avoid large vessels whenever possible.
3. Avoid the edge of the liver.

After selecting a detection location, mark with a bonding agent or wax crayon at the center position of the probe as the fibrosis detection location.







Fig. 7-4 Suitable localization

7.2.3.3 Image analysis of common errors

Image	Phenomenon	Causes and solutions
	1. The near field area is highlighted and blurred. 2. When moving the probe, you can see a blurry image of the fibrous membrane of the liver.	Probe position is too high, lung scan.

	1. Scanning of a bean-shaped object	Probe position is too low, kidney scan
	1. Several black shadows in the image.	Wrong probe direction, edge scanning.
	1. There are large vessels or tissues other than the liver, such as the bile ducts.	The position of the probe needs to be adjusted; you can move the probe along the intercostal space to avoid damage to non-liver tissue.

7.3 Elastography detection

7.3.1 Entering elastography mode

Press the "Mode" hot key to enter the elastography mode, as shown in Fig. 7-5.


Fig. 7-5 Elastography Mode Interface

7.3.2 Detection launch

After the ultrasonic coupling agent is applied to the probe, the clinician holds the probe to scan for fibrosis on the right side of the patient, as shown in Fig. 7-6. Hold the probe against the patient's skin with your right hand, with your right elbow pressed against the bed to prevent the probe from slipping, and hold the end of the probe with

your left hand to adjust the angle. (The distance between the operator, patient and other staff and the monitor should be $\geq 850\text{mm}$)

Fig. 7-6 Operation method

Note: Watch the pressure indicator panel in fig. 7-7 and keep the probe pressure in a moderate range.


Fig. 7-7 Pressure display panel

After the probe is in contact with the skin, as the probe pressure gradually increases, the pressure slider moves from left to right, and the color of the pressure indicator bar changes from pink to green and then to red, which means that the pressure changes from low to moderate pressure and then overpressure, respectively. The area between the green inverted triangle and the red inverted triangle on the pressure indicator bar is the moderate pressure area, and the blue box in this area is the highest-pressure area. To ensure accurate test results for patients with thicker subcutaneous fat, the pressure may

be slightly higher.

7.3.3 Mode M and mode A

The main part of the elastography mode interface is the black area in Fig. 7-5. It is noted that four blue rectangles are located in this area from left to right. The first two blue rectangles are used to display M mode and A mode. See fig. 7-8.


Fig. 7-8 M mode (left) and A mode (right) areas

Mode M and Mode A images can help determine the detected position. To obtain an accurate result, the liver fibrosis scanning algorithm requires that the liver tissue be relatively flat and at least 7 cm thick in the detected area. Fig. 7-9 shows an image of the M mode under three different conditions. (a) - image in M mode when the detected area is uneven. (b) - image in M mode when the detected area is uniform. And (c) - the image in mode M, when the detected area is uniform, and the examined patient holds his breath.


Fig. 7-9 Image in M mode under various conditions

1. Requirements for correct detection:

- 1) Mode A: the signal does not tremble strongly and quickly.
- 2) Mode A: There is no significant difference between different scan depths.
- 3) M mode: Image brightness of different depths is the same.
- 4) M mode: The image shows a normal stripe pattern.

Subcutaneous tissue, echo amplitude in mode A is large; M mode with high brightness, partially white cloth.

In the tissues of the liver, the amplitude of the echo signal in mode A is lower than in the subcutaneous tissues, and is relatively uniform, without large rises and falls at different depths. After the probe stabilizes, the Mode A signal becomes stable. In M mode, liver tissue appears dull over a wide range of depths, with no difference between depth and brightness. When the probe is stationary, tissues are presented in a layered form.

Tissues under the liver, the amplitude of the ultrasound echo signal in mode A and mode M are not similar to the liver tissue.

Fig. 7-10 Schematic diagram of correct detection

2. Error detection and root cause analysis

Schematic diagram					
Phenomenon	Mode M: There is no stable layer on the image	Mode M: Uneven image brightness at different depths	Mode M: Image gets darker and darker as depth increases	Mode M: The image becomes sharper and blurry as the depth increases	Mode M: There is no regular pattern in the image
Reasons	It is possible that the probe is not held stable. Probe may not be vertical to skin	On the path of the study, there are tissues that are not related to the liver.	The probe may touch a rib or the edge of a rib	Gas may be present in the path of inspection	Viewing position may not be correct

7.3.4 Elasticity modulus

Follow the M mode, identify a liver area with relative flatness and suitable thickness, ask the patient to hold their breath (to avoid the effect of breathing movements on the results) and press the footswitch. The results of this scan will be displayed on the interface. The shear wave velocity in units of m/s is shown in the lower right corner.




Fig. 7-11 Elasticity modulus

1. Correct detection

2. Error detection

Note:

1. It is better to use image scanning to find the best detection position before elastography detection, so as to improve the detection efficiency rate.
2. When the detection effect is unsatisfactory in the image scanning position (poor image quality or large relative deviation), slightly change the position and direction of the probe; if the adjustment is ineffective, switch to imaging mode to reselect the detection position.
3. When the effect of the breath on detection is high, you can ask the patient to hold their breath for 15 seconds each time and perform detection 3 times during each breath hold.
4. Detection is best done in the same position and at the same angle, as far as possible; please click the [Clear] button to clear the current records when you intend to change position or angle.
5. Reasons for a high-test result:
 - 1) The detection area is located at the edge of the liver.
 - 2) The edge of the probe is close to the ribs.

7.3.5 Detection results

When one study is completed with good data, the current number of detections and results will be updated in the right box. See Fig. 7-12.

Fig. 7-12 Successful detection

The internal algorithm will evaluate the inspection as unsuccessful if the position being checked is incorrect or the data is incorrect. The results will not be updated, and an error message will appear in the interface, as-shown in Fig. 7-13.

Fig. 7-13 Detection error

7.3.5.1 Stiffness

1. Stiffness is a parameter reflecting the level of liver fibrosis; the corresponding hardness value will be generated and displayed after each successful detection. Unit kPa. See Fig. 7-14.
2. The IQR reflects the degree of dispersion of the average 50% of the data. A lower value indicates a higher concentration; a larger value indicates a larger spread.
3. IQR/med indicates the percentage of absolute deviation from the mean for a single detection. IQR/med can only be used to weight the degree of deviation of a single detection result from the mean.

Fig. 7-14 Stiffness display section

Additional information: The detection result is invalid if IQR/med exceeds 33%.

7.3.5.2 Performance indicator

If the detection is valid, the system displays the Performance Score, Number of Valid Detections, and Total Detections. As shown in fig. 7-15 below, the total number of detections is 10, the number of valid detections is 10, so the displayed efficiency rate is 100%.

Fig. 7-15 Performance indicator display section

Additional information *If the success rate is less than 60%, this group of detection data is invalid and needs to be retested.*

7.3.5.3 Ultrasonic attenuation parameter

The unit of ultrasonic attenuation parameter (UAP) is dB/m. UAP may be used as an aid in the diagnosis and monitoring of patients with liver disease as part of a general liver assessment. The UAP is displayed as shown in fig. 7-16 while displaying IQR and IQR/med for UAP at the same time.

Fig. 7-16 UAP display section

Additional information: *The detection result is invalid if IQR/med exceeds 33%.*

7.3.5.4 Completion of detection

When the number of successful detections reaches 10 times, the survey ends. Click the [Up] and [Down] buttons to display all detection results, as shown in Fig. 7-17.

Measurements	
1. E =	14.1 kPa
2. E =	11.2 kPa
3. E =	10.2 kPa
4. E =	12.2 kPa
5. E =	10.4 kPa
6. E =	10.6 kPa
7. E =	10.4 kPa
8. E =	10.4 kPa
9. E =	10.4 kPa
10. E =	10.3 kPa

Fig. 7-17 Survey Results

Additional information: *The actual measurement results must comply with the following three factors:*

1. The number of measurements is at least 10 times.
2. The efficiency index exceeds 60%.
3. The IQR/med ratio is less than 33%.

7.4 Printing a report

Click the "Print" button and generate a liver fibrosis study report. The report contains information about the patient, the results of ten examinations and liver stiffness.

Chapter 8. Ultrasound Scan

8.1 Display modes

8.1.1 Mode B


Fig. 8-1 Mode B

To display in B mode only, press encoder [B]. In real time, it allows you to change image parameters by adjusting the menu bar. When paused, it allows you to switch on various measurement functions.

Function menu:

B Mode	Noise 0	Dyn 0	Enhancement 3	Chroma 0	Focus 9
	iClear High	Focus No. 1	Persist 0	Freq 3.0 MHz	Map 7
F1, F2	Noise, iClear High				
F3, F4	Dynamic, No. focus.				
F5, F6	Improvement, durability				
F7, F8	Colorize, Frequency				
F9, F10	Focus, Map				

Black and white reversion

Press [SHIFT+C] to return the black and white image to real time. Does not work in C mode.


Fig. 8-2 Black and white reversion

Up and down rotation

Press [SHIFT+ ()] or [SHIFT+>] to select the top or bottom image respectively in real time. Does not work in C mode.


Fig. 8-3 Up and down rotation

Rotation left and right

Press [SHIFT+<] or [SHIFT+>] to select the left or right image respectively in real time. Does not work in C mode.

Biopsy

Press [F7] on PC keyboard to enter real time biopsy.


Fig. 8-4 Biopsy

Scaling

In live mode, first turn the multi-dimensional adjustment knob to make the local zoom indicator light on, and then turn the multi-dimensional adjustment knob to adjust the zoom scale.

THI

Press [F8] on the PC keyboard to enter THI in real time. H3.5M will be displayed in the upper right corner of the screen.

Mark on the body

Press [F3] on the PC keyboard to display the mark on the body in real time.

Fig. 8-5 Mark on the body

X-ray cinematography

The Cine button will be displayed in the paused state.

Fig. 8-6 X-ray cinematography

8.1.2 Mode B/B


Fig. 8-7 B/B Mode

Press [Alt + M] to display the left image in B mode and freeze the right one; similarly, press [Alt + N] to display the right image in B mode and freeze the left. Changes image parameters in real time. Enables various measurement functions in the paused state. The function menu is the same as in B mode.

8.1.3 Mode 4B


Fig. 8-8 Mode 4B

Press [F1] on the PC keyboard to display images in 4B mode. Each click displays one image. Four clicks displays four images. The current image is just a live image. Changes image parameters in real time. Enables various measurement functions in the paused state.

The function menu is the same as in B mode.

8.1.4 Mode B+M


Fig. 8-9 Mode B+M

Press the [Alt + S] or [B/M] keys to display the image in B/M mode. In real time, this changes the parameters of the image; in pause mode, activates various measurement functions.

The function menu is the same as in B mode.

8.1.5 Mode M


Fig. 8-10 Mode M

Press the [Alt + R] or [B/M] keys in B+M mode to display the image in M mode. In real time, this changes the parameters of the image; in pause mode, activates various measurement functions.

The function menu is the same as in B mode.

8.1.6 Mode B+C


Fig. 8-11 Mode B+C

In B mode, press [Alt + H] or [CFM] key to turn on the B+C image display mode.

which changes image parameters in real time and enables various measurement functions in the paused state.

Function menu:

BC Mode	W.Filter 1	Ensize 13	PRF 3 kHz	Threshold 6	Baseline 3
	Speed Low	Freq 2.0 MHz	Noise 4	Persist 3	Map 0
F1, F2	W.Filter, Low speed				
F3, F4	Size, Frequency				
F5, F6	PRF, Noise				
F7, F8	Threshold, durability				
F9, F10	Baseline, Map				

8.1.7 Mode B+D


Fig. 8-12 Mode B+D

In B mode, press [Alt + G] or [PW] key to turn on the B+D image display mode, which changes image parameters in real time and enables various measurement functions in the paused state.

Function menu:

BD Mode	Mode Sync Syn	Gate 7	Volume 0	W.Filter 1	Scan Rate Mid
	Freq 2.0 MHz	Triplex sync Close	PRF 3 kHz	ChromaD 0	Baseline 3
F1, F2	Synchronization mode, frequency				
F3, F4	Gateway, Triplex Synchronization				
F5, F6	Volume, PRF				

F7, F8	W.Filter, ChrodmaD
F9, F10	Scanning speed, Baseline

8.1.8 Mode B+C+D


Fig. 8-13 Mode B+C+D.

Press [Alt + G] or [PW] in C mode to display B+C+D mode. It is possible to edit parameters in real time. In the pause state, the measuring function can be executed. The function menu is the same as in B+D mode.

8.2 General measurements and calculations

8.2.1 B, 2B General measurements

Range and accuracy of measurements of parameters: angle, volume, heart rate and time.

Name	Range of measurements	Accuracy of measurements
angle	$\pm 18^\circ$	$0^\circ \leq \text{angle} < 180^\circ$
Volume	$\pm 0.02874\text{m}^3$	$0 < \text{Volume} < 0.14137\text{m}^3$
Heart rate	$\pm 1 \text{ bpm}$	$40 \text{ bpm} < \text{Heart rate}$
Time	$\pm 0.575\text{s}$	$0\text{s} < \text{Time} < 11.5\text{s}$

Fig. 8-14 B/2B General measurements

8.2.1.1 Angle measurement


Fig. 8-15 Angle Measurement Diagram

Procedures:

1. In the pause state, press the [Measurement] key to open the measurement menu; Then use the trackball to move the cursor to "Angle" under "Measure C"; press the [enter] key and wait until the "+" cursor appears in the image area and the system enters the angle measurement state.
2. Use the trackball to move the cursor to point A; press the [Enter] key to commit.
3. Use the trackball to move the cursor to point O; press the [Enter] key to commit.
4. Use the trackball to move the cursor to point C; press the [Enter] key to commit.
5. Thus, calculate the included angle AOC (between segment OA and OC).
6. Repeat steps 1-4 to take another measurement.

Press the [Clear] key to clear the measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.1.2 Volume measurement

Clinically, the three-coordinate method is often used to determine the volume, the calculation formula is as follows:

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The calculation requires the measurement of three parameters (length, width and height) or area and distance based on the ultrasound image of the longitudinal and cross section. In B mode, first press [Alt + M] to display the left of the two images in B mode, and scan the patient to obtain a transverse (or longitudinal) section image; then press [Alt + N] to activate the right of the two images in B mode and scan the patient to obtain a longitudinal (or transverse) section image; freeze image.

This algorithm is applicable to measure the volume of the organs of an ellipse by measuring area and distance.

Procedures:

1. In the pause state, press the [Measurement] button to open the measurement menu; use the trackball to move the cursor to "Volume"; press the [Enter] button, the "+" cursor will appear in the imaging area and the system will enter the volume measurement state.
2. Use the trackball to move the cursor to one end of the ellipse's axis; press the [Enter] key to commit.
3. Use the trackball to move the cursor to the other end of the ellipse axis; press the [Enter] key to commit.
4. Use the trackball to change the length of the other axis.
5. Use the trackball to move the cursor to the start point of the distance; press the [Enter] key to commit.
6. Use the trackball to move the cursor to the end point of the distance; press the [Enter] key to commit.

Repeat steps 1-6 to take another measurement.

Press the [Clear] key to remove measurement marks

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.2 General measurement V+M, M


Fig. 8-16 General measurement B+M, M

8.2.2.1 Heart rate measurement

To improve the accuracy of measurements, this system requires the measurement of several heartbeat cycles and the calculation of their average value. The measurement start point means the starting points of two heartbeat cycles, and the measurement end point means the end points of two heartbeat cycles.

Procedures:

1. Use the trackball to move the cursor to the "HR" item; press the [Enter] key to enter the measurement state.
2. Use the trackball to move the cursor to the starting point and press the [Enter] key to commit.
3. Use the trackball to move the cursor to the end point and press the [Enter] key to commit.

Repeat steps 2 and 3 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.2.2 Valve speed measurement

This system calculates valve speed (tilt) by dividing the measured organ depth (distance) difference in B|M mode by the time interval. Different organ depths (distance) are expressed in millimeters; the time interval is expressed in seconds; the result is expressed in millimeters per second.

If the time interval approaches zero or the SL value is too high, the result will be displayed as "0".

Procedures:

1. In the pause state in B|M or M mode, press the [Measurement] key to open the measurement menu; use the trackball to move the cursor to "Vel"; press the [Enter] key to enter the measurement state.
2. Use the trackball to move the cursor to the starting point and press the [Enter] key to commit.
3. Use the trackball to move the cursor to the end point and press the [Enter] key to commit.

Repeat procedures 2 and 3 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3 General measurement PW

Fig. 8-17 General measurement PW

8.2.3.1 Time measurement

Measurement of the time interval between two points on a Doppler ultrasound image.

Procedures:

1. In PW pause state, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to "ET"; press the [Enter] key to enter the measurement state.
3. Use the trackball to move the cursor to the starting point and press the [Enter] key to commit.
4. Use the trackball to move the cursor to the end point and press the [Enter] key to commit.
5. After the calculation, the result will be displayed at the bottom of the imaging area. Repeat steps 2 and 4 to take another measurement.

Press the [Clear] key to delete the data and measurement symbol.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3.2 Heart rate measurement

Measure the time interval between two heartbeat cycles on a Doppler ultrasound image and calculate the heart rate.

Procedures:

1. In PW pause state, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to the "HR" item; press the [Enter] key to enter the measurement state.
3. Use the trackball to move the cursor to the starting point and press the [Enter] key to commit.
4. Use the trackball to move the cursor to the end point and press the [Enter] key to commit.
5. After the calculation, the result will be displayed in the imaging area. Repeat steps 2 and 4 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3.3 Speed measurement

Measurement of velocity and pressure drop at some point in the Doppler frequency spectrum on an ultrasound Doppler image.

Procedures:

1. In PW pause state, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to the "Vel" item; press the [Enter] key to enter the measurement state.
3. Use the trackball to move the cursor to the starting point and press the [Enter] key to commit.
4. Use the trackball to move the cursor to the end point and press the [Enter] key to commit.
5. After the calculation, the result will be displayed in the imaging area. Repeat steps 2-4 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3.4 Acceleration measurement

Measurement of velocities and time intervals at two points on a Doppler ultrasound image to calculate differential velocity and acceleration.

Procedures:

1. In the pause state and in PW mode, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to the "ACC" item; press the [Enter] key to enter the measurement state.
3. Use the trackball to move the cursor to the start points and press the [Enter] key to commit.
4. Use the trackball to move the cursor to the end points and press the [Enter] key to commit.
5. After the calculation, the result will be displayed at the bottom of the imaging area. Repeat steps 2-4 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3.5 General measurements

Dimension items:

PS: peak systolic velocity, namely the maximum velocity of erythrocytes in the sample volume.

ED: end diastolic velocity, namely end diastolic flow velocity.

TAMAX: The average flow rate of the maximum frequency spectrum.

RI: resistance index, $RI = RI = [(PS(m/c) - ED(m/c)) / PS(m/c)]$.

PPG: Peak differential pressure, pressure increment corresponding to peak systolic velocity, $PPG(mmHg) = 4 \times PS(m/s)^2$.

Procedures:

1. In PW pause state, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to the measured point and press the [Enter] key to fix PS.
3. Use the trackball to move the cursor to the measured point and press the [Enter] key to fix the MD.
4. Use the trackball to move the cursor to the measured point and press the [Enter] key to fix the ED.
5. After the calculation, the result will be displayed in the imaging area. Repeat steps 2-4 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.2.3.6 Doppler trace

Dimension items:

PS: peak systolic velocity, namely the maximum velocity of erythrocytes in the sample volume.

ED: end diastolic velocity, namely the flow rate at the end of the heartbeat cycles.

TAMAX: The average flow rate of the maximum frequency spectrum.

RI: resistance index, $RI = [(PS(m/c) - ED(m/c)) / PS(m/c)]$.

PPG: Peak differential pressure, pressure increment corresponding to peak systolic velocity, $PPG (mmHg) = 4 \times PS(m/s)^2$.

Procedures:

1. In PW pause state, press the [Measurement] key to open the measurement menu.
2. Use the trackball to move the cursor to the "Trace" item; press the [Enter] key to enter the measurement state.
3. Move the trackball to a path along the frequency wave.
4. Use the trackball to move the cursor to the starting points to be measured and press the [Enter] key to fix.
5. Use the trackball to move the cursor to the endpoints you want to measure and press the [Enter] key to lock.
6. After the calculation, the result will be displayed in the imaging area. Repeat steps 2-5 to take another measurement.

Press the [Clear] key to delete data and measurement marks.

During the measurement, the operator can press the [Pause] key to exit the measurement state.

8.3 Special measurements and calculations

8.3.1 Obstetric measurements

8.3.1.1 Measurement menu

Use the trackball to move the cursor to "Other" in the measurement menu. Select "Early OB Measurements" to enter the measurement menu.

8.3.1.2 Early OB Menu


Fig. 8-18 Early obstetric measurements

Select "Right Ovary" to enter the right ovary measurement menu.

Select "Left Ovary" to enter the left ovary measurement menu.

Select "Uterus" to enter the uterus measurement menu.

Select "Fetal Biometrics" to enter the fetal biometric measurement menu.

Fig. 8-19 Right ovary Fig. 8-20 Left ovary Fig. 8-21 Uterus Fig. 8-22 Fetal biometrics

8.3.1.3 OB measurement menu at later dates

Select "Fetal Biometrics" to enter the fetal biometric measurement menu.

Select "Longbones" to enter the Longbones measurement menu.

Select "AFI" to enter the amniotic fluid index measurement menu.

Fig. 8-23 Fetal Biometry. 8-24 Long bones Fig. 8-25 AFI Amniotic Fluid Index

8.3.2 Introduction to Obstetrics

Basic elements of measurement:

1. Choosing a default database version
2. Editing a user database
3. Calculation of gestational age and expected date of delivery
4. Fetal weight calculation
5. Calculation of gestational age and expected date of delivery based on last menstrual period
6. Calculation of gestational age and expected date of delivery based on basal body temperature
7. Assessment of fetal physiology
8. Fetal growth curve

8.3.2.1 Selecting a default database

The system provides four kinds of OB database versions for gestational age calculation:

- Asian version (Asia), European version (Europe), American version (America) total three versions of the fixed database;
- One editable version of the database (custom);
- For each obstetric version, parameters such as gestational sac (GS), biparietal head size (BPD), parieto-coccygeal length (CRL), femur length (FL), humerus length (HL), transverse abdomen size (TAD), vertebral length (LV), fronto-occipital dimension (OFD), abdominal circumference (AC), head circumference (HC) can be used to determine the duration of pregnancy and childbirth.

Form version Dimension element	Asia	Europe	America	Custom Version
GS	Tokyo	Hansman	Rempen	China
BPD	Tokyo	Headlock	Merz	China
CRL	Tokyo	Headlock	Rempen	China
FL	Tokyo	Headlock	Merz	China
HL	Osaka	Jinty	Jinty	Osaka
TAD	Hansman	Merz	Eriksen	Hansman
LV	Tokyo	Tokyo	Tokyo	Tokyo
OFD	Korea	Hansman	Chitty	Korea
Alternating current	Shinozuka	Headlock	Merz	Shinozuka
HC	Korea	Headlock	Merz	Korea

8.3.2.1.1 Editing a user database

The database can be edited by the user.

The editable database (UserCustom) uses statistical data applicable to the population of China. Don't make any adjustments without better statistics.

8.3.2.1.2 Calculation of gestational age and expected date of delivery

Weeks of pregnancy and date of delivery can be calculated from the obstetric sheet for measurement.

Parameters such as gestational sac (GS), biparietal head size (BPD), parieto-coccygeal length (CRL), femur length (FL), humerus length (HL), transverse abdominal size (TAD), vertebral length (LV), fronto-occipital dimension (OFD), collar space thickness (NT) are taken to measure distance, abdominal circumference (AC), head circumference (HC) is used by ellipse to measure perimeter | area; amniotic fluid (AFI) is used to measure four distances.

8.3.2.1.3 Assessment of fetal weight

User can measure AC, BPD, AC, FL or AC, FL, HC, then the data will be automatically stored in the internal memory, when calling the OB report, the system will automatically show the calculated result obtained by calculating several parameters in the report.

8.3.2.1.4 Calculation of gestational age and expected date of delivery based on last menstrual period

The general method for calculating due date was introduced in the previous chapter on automatic calculation of gestational age and due date using OB measurements. The following is a method for calculating weeks of pregnancy and due date based on last menstrual cycle data.

Regarding the method of manifestation during pregnancy, from the point of view of gynecology, the date of fertilization is taken as the first day, the entire gestational period is a total of 280 days, 10 months of pregnancy (28 days per month of pregnancy).

Since it is not easy for every pregnant woman to correctly estimate the time of pregnancy, for the clinical application of obstetrics, take the date of the beginning of the last menstrual period as the zero week calculation and zero date, this is the correct date of the period of delivery up to 40 weeks, zero day. In a pregnant woman, childbirth occurs within 38-42 weeks, this is a mature birth. This method of calculating is used for 28 days as a standard because menstrual weeks differ for each woman, so it is normal to have a 1-2 week discrepancy between the estimated due date and the actual due date.

<According to the last menstrual period calculation method>

$EDC(LMP) = LMP + 280(\text{days})$

(LMP: Date of last menstrual cycle, which means the first day of the last menstrual period)

(EDC: Estimated due date)

Procedures:

Move the trackball to select the "Basic Information" OB report page, click the drop-down button in the LMP (Last Menstrual Period) input field, select the date, then click the input field after "GA", it will display the estimated gestational age.

Fig. 8-26 Date selection dialog box

8.3.2.1.5 Calculation of gestational age and expected date of delivery based on basal body temperature

<Calculation method starting from ovulation date>

Estimated due date = ovulation date + 266(days)

Since it is difficult to determine the exact date of ovulation, usually the date of ovulation is determined by the date of transition from low temperature to high temperature at basal body temperature, i.e.

EDC (BBT) = BBT + 266 (days)

(BBT: Basal body temperature)

Procedures:

Move the trackball to select the OB report "Basic Information" page, click the drop-down button in the BBT (Basal Body Temperature) input field, select the date, then click the input field after "GA", it will display the estimated gestational age.

8.3.2.1.6 Biophysical assessment of the fetus

For biophysical assessment of the fetus, there are two types of common assessment methods:

Manning Scale Method — For each biophysical variation, score 2 for normal and 0 for exceptional. Observe the following parameters for 30 minutes: 1. Fetal breathing; 2. Fetal movement; 3. Muscular tone of the fetus; 4. Fetal heart rate; 5. Volume of amniotic fluid

The method of evaluation according to the Vincileos scale is the score "0", "1", "2" for each biophysical parameter.

1. Observe the fetal heart rate for 20 minutes; Observe the following parameters for 30 minutes: 2. Acceleration; 3. Fetal breathing; 4. Muscular tone of the fetus; 5. Volume of amniotic fluid; 6. Maturity of the placenta.

In clinical medicine, the Apgar scoring method is widely used to assess the condition of newborns, which includes the following five parameters:

1. **Skin coloration:** Assessment of the ratio between blood and oxygen for the lungs in newborns. 2 points - pink color of the whole body and extremities, 1 point - pink color of the body and cyanotic color of the extremities, 0 points - generalized pallor or generalized cyanosis.
2. **Heart rate:** Evaluation of the intensity and rhythm of the heartbeat in newborns. 2 points - good heartbeat > 100 times / min, 1 point - weak heartbeat < 100 times / min, 0 points - no heartbeat.
3. **Breath:** Evaluation of the maturity of the center and lungs in newborns. 2 points - normal breathing, 1 point - irregular breathing (for example: shallow, irregular or hurried and tense), 0 points - no breathing.
4. **Muscle tone:** Assessment of muscle tone in newborns. 2 points - active movements are expressed, 1 point - reduced muscle tone, some flexion of the limbs, 0 points - no muscle

tone, limbs hang down.

5. **Reflex excitability:** Evaluation of the ability to respond to external stimulation. 2 points - reaction in the form of movement, coughing, sneezing, loud cry, 1 point - the reaction is weakly expressed (grimace, movement), 0 points - no reaction.

Based on the above 5 parameters, a newborn is considered normal if he scores 10, a newborn who scores 7 or below needs closer observation, a newborn who scores 4 or below needs closer observation, most newborns score between 7-10, the doctor conducts appropriate treatment according to the child's score. After clearing the airways and inhaling oxygen, the newborn soon gets better.

Clinical values: A high or low fetal score is highly correlated with neonatal Apgar scores, fetal anxiety, and antenatal and perinatal mortality. The lower the score, the higher the likelihood of such factors. When all assessed parameters are normal, "10" points corresponds to fetal physiology, perinatal mortality is 0; conversely, if all parameters are abnormal, the fetal physiology score is "0" points, fetal mortality is 40%.

The Winzelos score is used in the system to get a score for the following parameters:

- Fetal movements
- Muscular tone of the fetus;
- Fetal breathing
- Volume of amniotic fluid;
- No voltage test;
- Evaluation of the placenta.

Procedures:

Move the trackball to select the OB report "Calculation" page, click the drop-down button in the "Fetal Movements" input field, select the desired score, then repeat the operation for other parameters.

8.3.2.1.7 Fetal growth curve

The fetal growth curve is displayed on the growth curve page of the OB report.

After entering the OB Trend page, click the ten items GS, BPD, CRL, FL, HL, TAD, LV, OFD, AC, and HC, which refer to the ten corresponding growth curves.

Select "Picture Option" to display the menu for adjusting picture settings:


Fig. 8-29 Set the default value of parameter 1


Fig. 8-30 Set the default value of parameter 2

Use the trackball and press the [Enter] button to set the parameters to the default value.

8.5 Printing a report

8.5.1 Saving/accessing an image

Save the ultrasound image if necessary. To do this, click the "Snapshot" button on the ultrasound interface to automatically save the current ultrasound image, as shown in

Fig. 8-31.

Fig. 8-31 Saving an ultrasound image

The saved image will appear at the bottom of the screen. It is possible to save multiple consecutive images as shown in fig. 8-32.


Fig. 8-32 Saving multiple ultrasound images

If you want to access any of the saved images, double-click on such an image. The image will appear in the middle of the screen as shown in Fig. 8-33.

Figure 8-33 Viewing the selected ultrasound image by double-clicking

8.5.2 Printing an image scan report

Printing a diagnostic report of the general examination in A4 format on an HP color inkjet printer. The user can check the "Report" box and select an image to be added to the diagnostic report. The report should contain 1-4 ultrasound images and a place for the doctor to enter "Comment", as shown in fig. 8-34.

In the ultrasound imaging scan interface, click the "Print" button to print the ultrasound report. The next paragraph will describe the printing process.

1. First, check that there is paper in the printer and that the printer is working properly. The printout is only available if the printer is working.

- 1) Good condition indication. The green light on the power switch remains on.
- 2) A sign of a malfunction. The printer's ink status light and the Recovery light are on, or the Power light is flashing.
- 3) In this case, correct the problems before printing. See the HP Printer User Guide for more information.

2. Click the "Print" button in the ultrasound imaging interface. Wait for the pop-up dialog box to appear.

3. The printer will start printing a diagnostic report. The user can perform other operations after closing the dialog box.

4. The printer exports the printed report. By then, the printing will be complete. The process takes about 2 minutes.

⚠ Warning: Do not connect or disconnect the power cable and signal cable while the printer is running. Otherwise, it may cause damage or malfunction of the system and printer.


Fig. 8-34 Ultrasound Image Report

8.5.3 Saving a report

Click the "Save"  button on the report interface to save the report.

Select the export path in the pop-up dialog box, as shown in Fig. 8-35, enter a file name and click the "Save" button to finish saving the report. The report format is JPG.

⚠ Hint: You can save the report to a local drive, an external USB drive, or an external mobile hard drive.


Fig. 8-35 Save the report


8.6 Working with shared buttons


Click on the button and navigate to the parameter to be configured. The indicator shows that the function is selected. You can adjust settings such as Depth, Focus and Zoom. Turn the knob to adjust the corresponding parameters, clockwise to increase, and counterclockwise to decrease the parameters.

This button allows you to scroll and select menus.

Click this button to quickly enter imaging mode. It can provide gain adjustment, the eight-stage gain potentiometer can also achieve the function of gain adjustment.

Scroll through the images of the probe you are using.

[Cancel]

Press the [Cancel] key on the control panel or the Cancel button on the interface to cancel the current detection. If there is a detection record, the system will automatically save the result and exit to the detection interface. Cancel button as shown in Figure 8-36.

Fig. 8-36 Cancel Button

[Snapshot]

Press the [Snapshot] button on the control panel to take snapshots of ultrasound images in real time. Snapshot effects as shown in fig. 8-37.

Fig. 8-37 Snapshot effects

[Save]

Press the "Save" button on the control panel or the "Save" button on the screen to save the diagnostic information and exit. The "Save" button, as shown in fig. 8-38.

Fig. 8-38 Save Button

[Up]/

[Down]

Press the [Up] and [Down] buttons or the "Up" and "Down" buttons on the interface to view the real-time ultrasound images. The "Up" and "Down" buttons are shown in fig. 8-39.


Fig. 8-39 Up and Down Buttons

[Print]

Press the [Print] button on the control panel or the "Print" button on the interface to print the saved images. The "Print" button, as shown in fig. 8-40.



Fig. 8-40 Print Button

[Pause]

Press the [Pause] button on the control panel to pause/resume the real-time ultrasound image, and then the pause mark will be displayed on the interface, as shown in Fig. 8-41.

Fig. 8-41 Pause mark

Single mode

Check the "Single" box as shown in fig. 8-42 to set single mode and print the last selected ultrasound image. If you do not select single mode, the system will automatically print the two selected ultrasound images.


Fig. 8-42 "Single" option

Freezing the image

The "Freezing the image" button on the interface performs different functions in different states. The "Freezing the image" button, as shown in Figure 8-43.

- 1) When the ultrasound image is not a snapshot, the Freezing the image function also performs the Snapshot function.
- 2) If the ultrasound image was captured using the Snapshot button, click Freezing the image to stop scanning the image and preview the image in real-time magnification.


Fig. 8-43 Freezing the image button.

Image launch

In the scan stop state, click "Image launch" on the screen, and the system will re-enter the image scan state. "Image launch" button, as shown in Figure 8-44.


Fig. 8-44 Image launch button.

[Clear]


Press the [Clear] button on the control panel to clear the ultrasound image measurement data and image measurement lines.

Cleaning

Click the "Clear" button to restart the image scan. Clear button as shown in Figure 8-45.

Fig. 8-45 Clear Button

[Measurement]

Adjust the image, press the [Pause] button on the control panel to unfreeze the image, and then press the [Measure] button on the control panel to enter the measurement mode.

Chapter 9. Work after use

9.1 Turning off the main unit

1. Turn off the main unit on the control panel or use the "Power off" button on the control interface after use.
2. Turn off the monitor.
3. Turn off the main power.

9.2 Power off

You must turn off the power source after each use. Please check the following points in advance:

1. The probes are placed in an appropriate holder, and the cable is on the stand.
2. The power of the connected peripheral equipment has been turned off in accordance with the instructions in the manual.

The system power switch is located on the rear of the machine.

9.3 Checking of the technical condition

After each use, a preventive inspection should be carried out. The check points and procedure are the same as the check before use.

⚠ Warning: *If any deviations are found, please contact our sales office, customer service department or agent.*

Chapter 10. Acoustic power

 **Warning:** *This section of the main user guide applies to the entire system, including the main unit, probes, accessories, and peripherals. This section contains important safety information for product operators regarding acoustic output and how to control patient exposure using the ALARA (Least Possible Impact Achievement) principle. This section also contains information about acoustic power testing and real-time output power display.*

The system does not create noise during operation.

Please read this information carefully before using the system.

10.1 Problems related to biological effects

Diagnostic ultrasound is considered safe. In fact, there have been no reports of patient injuries caused by diagnostic ultrasound.

It cannot be exactly stated that ultrasound is 100% safe. Studies have shown that extremely high intensity ultrasound is harmful to body tissues.

Over the past few years, diagnostic ultrasound technology has made great strides forward. This rapid progress has raised concerns about the potential risk of biological impact when new applications or diagnostic technologies emerge.

10.2 Reasonable use statement

Although there are no confirmed biological effects on patients caused by exposure to existing diagnostic ultrasound instruments, there is a possibility that such biological effects may be identified in the future. Therefore, ultrasound should be used with caution to ensure a medical benefit to the patient. High exposure levels and long exposure times should be avoided when appropriate clinical information is obtained.

10.3 ALARA principle (achieving the result with the lowest possible impact)

When using ultrasonic energy, the ALARA principle must be applied. The application of the ALARA principle ensures that the overall energy level is controlled below a level low enough that no biological effects are generated during the accumulation of diagnostic information. The total energy is controlled by the output power and the total emission time. The intensity of emission required for the examination varies depending on the patient and the clinical case.

Not all examinations can be performed with extremely low levels of acoustic energy.

Controlling the acoustic level at an extremely low level results in poor image quality or an insufficient number of Doppler signals, which negatively affects the reliability of the diagnosis. However, an increase in acoustic power above the required value does not always improve the quality of information needed for diagnosis, but rather increases the risk of biological effects.

Users must be responsible for patient safety and use ultrasound conscientiously. Conscious use of ultrasound means that the output power of the ultrasound must be selected based on ALARA.

Additional information on the concept of ALARA and the possible biological effects of ultrasound is available in the AIUM (American Institute of Ultrasound Medicine) document titled "Medical Ultrasound Safety".

10.4 Explanation of MI/TI

10.4.1 Basic information about MI and TI

The relationship of various output parameters of ultrasound (frequency, acoustic pressure and intensity, etc.) with biological effects is currently not fully understood. It is recognized that two main mechanisms can cause biological effects. One is a thermobiological effect with the absorption of ultrasound by tissues, and the other is a mechanical biological effect based on cavitation. The Thermal Index (TI) refers to the relative index of temperature rise due to thermobiological forcing, and the Mechanical Index (MI) refers to the relative index of mechanical biological forcing. The TI and MI indexes reflect instantaneous output conditions, so they do NOT take into account the cumulative effects of total test time. The TI and MI models contain practical simplifications for the interaction of complex biological effects. The operator must be aware that the actual temperature rise may, in the worst case, be several times higher than the displayed TI value.

MI (Mechanical Index):

Mechanical biological effects are the result of compression and decompression of isolated tissues with the formation of microbubbles, which are called cavitations.

MI is an index that indicates the possibility of cavitation formation based on acoustic pressure, and the value at which the peak acoustic vacuum pressure is divided by the square root of the frequency. Therefore, the MI value becomes smaller when the frequency is raised or the peak rarefied acoustic pressure is lowered, and cavitation becomes difficult to occur.

$$MI = \frac{Pr, \alpha}{\sqrt{f_{\text{ant}}} \times CMI}$$

$$CMI = 1 \text{ (MPa}/\sqrt{\text{MHz}})$$

For a frequency of 1 MHz and a peak vacuum acoustic pressure of 1 MPa, MI becomes equal to 1. It can be considered that MI is one of the cavitation generation thresholds. It is especially important to maintain a low MI value when both gases and soft tissues co-exist, for example when treating the lungs on a heart scan and intestinal gases on an abdominal scan.

TI (Thermal index):

TI is defined as the ratio of total acoustic power to the acoustic power required to raise tissue temperature by 1 °C. In addition, since the temperature increase varies greatly depending on the tissue structure, TI is divided into three types: TIS (Soft Tissue Thermal Index), TIB (Bone Thermal Index) and TIC (Skull Thermal Index).

TIS: Thermal index referring to soft tissues in areas such as the abdomen and heart area.

TIB: Thermal index for use in areas such as the fetal (second and third trimester) or neonatal brain (via the fontanel), where the ultrasound beam passes through soft tissue and the focal area is in close proximity to bone.

TIC: Thermal index for use in areas such as pediatric and adult skulls, in which the ultrasonic beam passes through the bone near the entrance of the beam into the body.

Although power output is automatically adjusted for selected applications, high TI values should be minimized or avoided in obstetrics. The WFUMB (World Federation of Ultrasound in Medicine and Biology) guidelines state that a temperature increase of 4 degrees Celsius for 5 minutes or more should be considered potentially damaging to fetal tissue.

The lower the MI/TI values, the lower the biological effect.

10.4.2 MI/TI display

The TI and MI values are displayed in real time at the top of the screen. The operator must monitor these index values during examinations and ensure that exposure times and output values are kept to the minimum required for effective diagnosis. The display accuracy is 0.1.

Under various operating conditions, as soon as the situation occurs when the MI value exceeds 1.0, the starting point of the display of MI values is 0.4.

In the same way, once a situation occurs where the TI value is greater than 1.0, the TI value will be displayed and the starting point is 0.4. You can specify which TI element should be displayed in the settings.

Note: *If the MI or TI value is greater than 1.0, you must be careful to follow the ALARA principle.*

10.5 Acoustic power management

A skilled operator can use the system controls to limit ultrasound power and adjust image quality. There are three categories of system controls in terms of power. These are controls that directly affect power, controls that control power indirectly, and controls that are receiver controls.

Direct control:

If necessary, you can control the acoustic power using the "Quick Menu" control on the control panel. In this case, the maximum acoustic power never exceeds MI 1.9, TI 6 and ISPTA.3 720 MW/cm² in any mode of operation.

Indirect control:

Controls that indirectly affect power are many image parameters. These are operating modes, frequency, focal point position, image depth and pulse repetition rate (PRF).

The mode of operation determines whether the ultrasonic beam is scanning or non-scanning. The thermobiological effect is closely related to the M mode, Doppler and color modes.

The acoustic damping of tissue is directly related to the frequency of the transducer.

The focal point is related to the active aperture of the transducer and the beamwidth.

The higher the PRF (pulse repetition frequency), the more output pulses occur in a certain period of time.

Receiver controls:

Receiver controls (e.g. gain, dynamic range, image post-processing, etc.) do not affect power. These controls should be used whenever possible to improve image quality before using controls that directly or indirectly affect power.

10.6 Acoustic power limits

As required by FDA Track 3, a reduction approach has been included within the FDA acoustic power limits as noted below. The maximum acoustic power level of any transducer in any operating mode is expected to be below these limits.

FDA Maximum Sound Power Limits for Track 3 (Reduced Values)

Application	$I_{spta.3}$ (MW/cm ²)	$I_{sppa.3}$ (W/cm ²)	MI
Areas ☐ (except eyes)	720	≤ 190	≤ 1.9

Chapter 11. System Maintenance

11.1 Probe Maintenance

⚠ Warning: *Do not service the probe while a patient is using it.*

11.1.1 Precautionary measures

Degree of protection of the probe for scanning fibrosis: IPX4.

Degree of protection of the image probe: IPX4.

Note:

1. *Probe protection. Do not strike the probes or let them fall.*
2. *Use a recognized ultrasound gel.*
3. *Do not connect or disconnect probes during operation.*
4. *Do not bend or pull on the probe cable, power cable, or other items.*
5. *Do not allow probes to come into contact with paint thinner, ethylene oxide, or other organic solvents.*
6. *The cleaning method and procedure must comply with the instructions.*
7. *Periodically check the probe cable, socket, and acoustic window area.*
8. *Turn the power off and off before connecting or disconnecting the probe.*
9. *Do not drop the probe on the floor or on a hard object, avoid colliding with the acoustic window of the transducer to avoid damaging it.*
10. *Place the probe in the probe holder when not in use.*
11. *Do not expose the probe to external heat.*
12. *Do not bend or stretch the probe cable, otherwise the internal connection may break.*
13. *Only use bonding agent on the probe head, clean the probe after use.*
14. *The Fibrosis Scanning Probe is recommended to be serviced, calibrated, or replaced after repeated use to avoid measurement drift.*
15. *Do not turn the probe sonic head to scan for fibrosis.*
16. *Only press the sonic head of the fibrosis probe during the examination. Inside the probe there are pressure-sensitive components. You may damage the probe if you apply too much pressure.*

11.1.2 Storage

1. *Place the probes in a clean and dry environment away from direct sunlight.*
2. *Store probes at -5°C to +40°C. Do not place probes in high pressure or vacuum environments.*
3. *Handle probes with care and avoid damage during storage.*

4. Store the probes in the probe case during transport or when not in use.

11.1.3 Cleaning and disinfection

⚠ Warning: Do not immerse the probe connector in liquid such as water or disinfectant. Immersion may result in electric shock or malfunction.

⚠ Attention!

1. When cleaning the device and probes, be careful not to damage the important marks showing information such as "Machine Model", "Machine Serial Number", "Monitor Model and Serial Number", "Probe Serial Number" and "Manufacturing Date" on probes.
2. Wear sterile gloves when cleaning and disinfecting the transducer to prevent infection.
3. After disinfection, rinse the probe thoroughly with sterile water to remove any chemical residue. Chemical residues on the probe can be harmful to the human body.
4. HISKY does not guarantee the effectiveness of the disinfectant and sterilization solution. Contact the respective manufacturers for more details.
5. If you are using the supplied solution, pay attention to the period of use of the solution, make sure that it has not expired.
6. Determine the degree of disinfection depending on the type of tissue in contact during use.
7. According to the clinical application of the probe to determine the appropriate concentration and cleaning time. For more information, see the instructions on the disinfectant label.
8. The fibrosis scanning probe is not sealed, do not use water, soaking, steam or other cleaning and disinfection methods, which may cause the disinfectant to leak into the internal system and damage the probe.
9. Do not use hydrocarbon-based glass cleaners to clean the monitor.

Note:

1. Wipe off the ultrasound gel thoroughly after the examination. Otherwise, the ultrasound gel may harden and degrade the image quality of the probe.
2. DO NOT allow the probe to overheat (more than 55°C) during cleaning and disinfection. High temperatures can deform or damage the probe.
3. Do not use autoclave, ultraviolet rays, gamma rays, gas, steam or heat sterilization methods to disinfect the transducer, otherwise serious damage may result.
4. Do not use iodine, mercury or organic mercury chloride, organic solvents, peracetic acid, acid solution, alkaline solution to clean the probe.
5. The use of cationic surfactants such as benzalkonium bromide, acetate,

chlorhexidine, poly(hexamethylene biguanide) hydrochloride is not recommended. It is not recommended to use nano silver inorganic antimicrobial agents and ethylene glycol with toxicity, so as not to affect the ultrasonic transmission. Do not use disinfectant gel containing cellulose.

11.1.3.1 Cleaning and disinfecting the imaging Probe

1. Wear sterile gloves to prevent infection.
2. Remove any binding gel from the probe by wiping it with a sterile cloth or gauze.
3. Wipe surfaces with 75% rubbing alcohol.
4. Dry the probe with a sterile cloth or gauze. Do not dry the probe by heating it.

11.1.3.2 Cleaning the Fibrosis Scan Probe

1. Wear sterile gloves to prevent infection.
2. Remove any binding gel from the probe by wiping it with a sterile cloth or gauze.
3. Wipe surfaces with 75% rubbing alcohol.
4. Dry the probe with a sterile cloth or gauze. Do not dry the probe by heating it.

11.1.3.3 Cleaning the trackball

1. Turn off the system and unplug the power cord from the outlet.
2. Place your fingers on the guide point of the retaining ring.
3. Slowly turn the retaining ring counterclockwise.
4. Remove the retaining ring and guide ball from the keyboard.
5. Wipe off oil and dust from the guide ball with a dry cloth.
6. Remove oil and dust from the retaining ring with a cotton swab.

11.1.3.4 Cleaning of the control panel

1. Wipe the surface of the control panel with a lint-free cloth lightly dampened with a mild cleaner.
2. Wipe the buttons and keyboard with a cotton swab.
3. Remove any solids between the buttons with a toothpick.
4. Remove moisture with a dry soft cloth and allow hard surfaces to air dry completely.

Note: *Do not spill water or other liquid on the system during cleaning.*

11.1.3.5 Cleaning of the footswitch

1. Wipe the surface of the footswitch with a lint-free cloth lightly dampened with a mild cleaner.
2. After cleaning, wipe the surface with a clean, lint-free cloth.

11.1.3.6 Cable cleaning

1. Wipe the surface of the cables with a lint-free cloth lightly dampened with a mild detergent.
2. After cleaning, wipe the surface with a clean, lint-free cloth.

11.1.3.7 Cleaning of the main unit

1. Turn off the system.
2. Turn off the power.
3. Wipe the surface of the main unit with a lint-free cloth. Or you can use a lint-free cloth lightly dampened with a mild detergent.
4. After cleaning, wipe the surface with a clean, lint-free cloth.

Note: Do not spill water or other liquid on the main unit, the main unit is IPX0.

11.2 Maintenance of the main unit

Warning:

1. Before cleaning the system, turn off the power source and unplug the power cord from the outlet. Cleaning a live system creates a risk of electric shock.
2. Do not perform maintenance on the main unit while the patient is using it.

11.1.2 Precautionary measures

1. Ensure good working conditions and keep equipment clean.
2. Restart the main unit at least 5 seconds after stopping.
3. Protect the control panel film. Use the control panel carefully and avoid scratching or tearing the film with sharp objects.
4. Wipe the main unit with a soft cloth dampened with water. If it is still dirty, clean the main unit with a damp soft cloth and detergent. Do not splash water into the main unit.
5. Do not clean the monitor with a hydrocarbon-based glass cleaner.

11.2.2 Maintenance

1. When using and cleaning the system, try to protect as much as possible important marks such as "Machine model", "Main unit number", "Probe number", and "Delivery date".
2. Keep the user manual carefully if it is not used frequently.
3. Check the system periodically. The following paragraph describes the detailed method.

11.2.3 Daily checklist

1. After each use, completely clean the gel on the probe with a soft cloth dampened with water.
2. Make sure the probes and cable are not damaged or aged. Stop using any such probes and cables immediately.
3. Clean the used probe.

⚠ Warning: *Stop using any probe, cable, or other damaged item immediately. Otherwise, it will lead to the danger of electric shock.*

11.2.4 Weekly checklist

1. Check probes for damage.
2. Make sure the power cable is not frayed or damaged.
3. Disinfect the used probe.
4. Turn off the ultrasonic power supply B, remove the plug from the socket and clean.

⚠ Hint

1. *Turn on the system and check it every three months if it has been turned off for a long time.*
2. *The battery of the main board should be removed by our technician, not by the user, if the system is in a power off state for a long time.*

11.3 Security check

11.3.1 Security check

Qualified personnel should carry out the following safety checks at least once every 24 months, contact your local service representative for periodic maintenance.

1. Check the system and accessories for mechanical and functional damage.
2. Check if the safety labels can be easily identified.
3. Make sure the fuse is suitable for rated current and no-load characteristics.
4. Check the protective earth resistance. The impedance value must not exceed 0.1 ohm.
5. Check earth leakage current: $\leq 500 \mu\text{A}$ under normal conditions;
 $\leq 1000 \mu\text{A}$ in case of one failure.
6. Check leakage current per case: $\leq 100 \mu\text{A}$ under normal conditions;
 $\leq 500 \mu\text{A}$ in case of one failure.
7. Check patient leakage current:
 $\leq 100 \mu\text{A}$ under normal AC conditions;

- ≤ 500 μ A in the event of a single AC failure;
- ≤ 10 μ A under normal DC conditions;
- ≤ 50 μ A in the event of a single DC failure;
- ≤ 100 μ A under normal conditions;
- ≤ 500 μ A in case of one failure.

8. Check patient leakage current (apply 110% supply voltage to the working part to check its operation in case of one failure): ≤ 5 mA.

If any of the security checks fail, stop using the system immediately. Contact the manufacturer or an authorized service provider for repairs.

11.3.2 Troubleshooting

No	Fault	Handling
1	The power indicator lamp does not light after the power switch is turned on.	Check the power supply system, fuse and cable.
2	The power indicator light is on, but the monitor does not show anything.	1. Reboot the system 5 minutes after stopping. 2. Adjust the monitor's contrast and brightness controls. 3. Check or re-insert the data row.
3	Ultrasound A or Ultrasound M is not shown on the display.	Make sure the probes are connected correctly.
4	The image is not displayed correctly and contains noise.	1. Remove electrical equipment. 2. Do not turn on electrical equipment when turning on the ultrasound system.
5	The probe is not shown on the display.	1. Press the [Probe] button to check if the system recognizes the probe; otherwise, replace the probe socket.
6	Image interference (longitudinal deformation).	1. Check if the earth is connected to the socket. 2. Power off the printer. 3. Find the source of network interference (starting or stopping nearby large equipment).
7	There is no ultrasound output in M and A modes.	1. To resume testing, press the [Pause] button twice. 2. Press "Mode" button to switch to imaging mode, press "Mode" button again to return to fibrosis scanning mode. 3. In the paused state, reconnect the probe to scan for fibrosis.
8	Footswitch does not work or does not respond.	1. Make sure the probe pressure is in the temperate zone. 2. Check if the probe tip is lubricated with a coupling agent. 3. Reconnect the footswitch.

9	The presence of pressure when there is no load on the probe (the pressure indicator lights up pink or green).	1. Lightly pull the black tab on the probe.
10	The trackball doesn't work well.	1. Connect a wired mouse, if available. 2. Open the trackball cover, remove the guide ball, clean the internal dust, put the guide ball back in place.
11	The keyboard is not available or input is not working correctly.	1. Make sure the three little blue lights are on in the top right corner. 2. Press the Caps Lock key on your keyboard and see if the indicator changes.
12	Action buttons do not work correctly.	1. Check if a single button is disabled or if all buttons are disabled. 2. Press the button on the screen to check if it works correctly. 3. Check if it is in the correct mode of operation.
13	Lack of response of buttons, incorrect display of the image and other functions when the above problems cannot be solved.	1. Press the probe button, reset the probe. 2. Press the [Pause] button several times to check if the image can be switched between freeze and defrost. 3. In the scanning pause state, press the "Clear" button on the screen, wait three seconds to restart the imaging module, do not press again. 4. Attempt 3 is invalid, click "Logout" and try to log in again. 5. Attempt 4 is invalid, shut down the system and restart after 5 minutes.

⚠ Attention! *The manufacturer can provide relevant wiring diagrams and other technical information according to customers' requirements.*

⚠ Warning:

1. *Make sure the system is turned off before changing the fuse. When replacing a fuse, the correct type must be selected to avoid mismatch. If the system still does not work, turn off the power, contact the manufacturer or an authorized service center. Do not open the machine for repairs without permission to prevent personal injury and property damage.*
2. *If our special filtered power cable is used, be sure to unplug it from the outlet in time after completing the test and turning off the product.*
3. *Our special power cable or filtered socket can only be connected to our machine. Non-intended use is prohibited.*
4. *In case the failure of the machine cannot be resolved with the above measures, contact the Manufacturer or an authorized service center, providing the necessary information: description of the failure, machine serial number, probe serial number,*

date of purchase.

5. *Two spare fuses are supplied with the machine, make sure the main power supply is switched off when replacing.*
6. *If a spare fuse is used, the user is not allowed to purchase a fuse other than the T6.3AH250V fuse.*
7. *If replacing the fuse does not resolve the service problem, turn off the main power immediately and contact the manufacturer or an authorized service center. Do not, under any circumstances, open the cover of the machine without permission to avoid personal injury and/or property damage.*

11.4 System repair

In the event of a system malfunction, turn off the power supply immediately. Tell our specialists the detailed information about the problem so that they can provide the service in time. Detailed information about anomalies will greatly facilitate our maintenance work and save maintenance time.

If the fuse is blown, replace it with a T6.3AH250V fuse according to the instructions of our technicians. The scope of delivery includes two fuses of the specified specifications. Such fuses are also available on the market. If the fuse blows again after replacement, notify our technicians immediately.

This system contains consumables and parts that need to be replaced, such as imaging transducer, fibrosis scanning transducer, footswitch, power cable, etc. Before replacement, contact the manufacturer for instructions.

The system is a complex electronic system. If the system needs to be repaired or adjusted, please contact us. We will send a service technician to repair the system at your site. We are not responsible for any malfunctions or hazards associated with the repair of the system by others. Hisky will send professional service personnel trained and authorized to perform maintenance.

11.5 Disposal and environmental protection

Do not dispose of waste, residues and other substances generated during the use of the system. Follow the rules for the disposal of medical waste. Do not use or disassemble the system or accessories at the end of their service life.

When decommissioning the product, the user must disinfect it (according to current environmental regulations, an undisinfected product is classified as hazardous waste).

There are three ways to dispose of the product:

1. Disposal by the manufacturer.
2. Disposal by a company licensed to collect equipment for the purpose of destruction or disposal, subject to all necessary measures to protect the life and health of people and the environment.

3. Recycling by the user, if he has the appropriate technology and resources.

Do not dispose of this product in the trash.

Send a recycling notice to the manufacturer.

11.6. Packing

There is no consumer packaging for the device.

The transport packaging of the device must comply with the requirements of regulatory documentation and protect the product from contamination, mechanical damage, climatic influences during transportation, handling and storage.

The monitor must be packed in a cardboard box with foam inserts and the cardboard box must be sealed with packing tape.

Probes for scanning and imaging should be packed in polyethylene bags and packed in a separate case.

All components of the product must be packed in polyethylene bags and then, together with a case with probes, they are placed in a cardboard box with foam nests and sealed with packing adhesive tape.

The main unit should be packed in a waterproof nylon bag and then packed in a PE bag.

The instruction manual must be enclosed in a plastic bag.

Then, the monitor in a cardboard box, all components together with the case, packed in a cardboard box, the main unit and the instruction manual are placed in a transport container - a wooden box.

Overall size of the transport package, see $L \times W \times H = 110 \text{ cm} \times 75 \text{ cm} \times 135 \text{ cm}$.

As packaging elements, use packaging materials - according to the regulatory documentation for a specific material, providing the necessary conditions for transportation and storage.

The gross weight of one packing piece must be: $130 \text{ kg} \pm 10\%$.

The marking must be clearly legible and retained throughout the entire period of transportation and storage.

11.7. Data on the labeling packaging of a medical device

Markings and inscriptions must be applied to each medical product "Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis" supplied to the consumer.

Marking must be applied to the labels in accordance with international standards and Russian standards: IEC60601-1:2005/A1:2012, IEC60601-1-2:2014, ISO 15223-1:2016, ISO 780: 2015, 15223-1:2020.

On the front side of the body of the Device, a plate must be attached, on which the trademark of the manufacturing organization "iLivTouch" must be indicated.

Note: Hereinafter, "manufacturer's trademark" is a trademark of "iLivTouch" (International Registration No. 1501683 (ILIVTOUCH)).

On the plate on the back side of the body of the device and on the transport label, the trade organization-manufacturer "HISKY MED" must be indicated.

11.7.1 Main unit marking

11.7.1.1 Manufacturer's mark, model: Pro1000.

Shear wave quantification ultrasound diagnostic system

Model: Pro1000

Power: 100-240V~ 50/60Hz 350V-A Weight: 81kr

Company name: Wuxi Hisky Medical Technologies Co., Ltd.

Address: Office B401, 530 Plaza, University Science Park, Taihu International Science and Technology Park, 214135 Wuxi, Jiangsu Province, PEOPLE'S REPUBLIC OF CHINA

European representative

Company name: Shanghai International Holding Corp. GmbH (Europe).

Address: Eiffestrasse 80, 20537 Hamburg, Germany.

Made in China.

Fig. 11-1 - Marking of the plate of the main unit model Pro1000 on the back of the machine.

11.7.1.2 Manufacturer's mark, model: Pro3000 advanced search

Shear wave quantification ultrasound diagnostic system

Model: Pro3000 advanced search

Power: 100-240V~ 50/60Hz 350V-A Weight: 81kr

Company name: Wuxi Hisky Medical Technologies Co., Ltd.

Address: Office B401, 530 Plaza, University Science Park, Taihu International Science and Technology Park, 214135 Wuxi, Jiangsu Province, PEOPLE'S REPUBLIC OF CHINA

European representative.

Company name: Shanghai International Holding Corp. GmbH (Europe).

Address: Eiffestrasse 80, 20537 Hamburg, Germany.

Made in China.

Fig. 11-2 - Marking of the plate of the main unit model Pro3000 advanced search on the back of the machine.

11.7.2 Marking for delivery to the Russian Federation

A plate must be affixed to the back of the apparatus housing, indicating:

- name and model of the medical device;
- full and abbreviated (if any) name of the manufacturer;
- manufacturer's postal address;
- full and abbreviated (if any) name of the authorized representative of the manufacturer;
- postal address of an authorized representative of the manufacturer;
- trademark of the manufacturer;
- "Serial number" symbol and related data;
- symbol date of manufacture and corresponding data (month, year);
- nominal supply voltage range;
- rated power frequency;

- Power consumption;
- class I product protection symbol;
- CE conformity mark;
- symbol "Refer to the instructions for use";
- instructions for disposal in the form of a symbol;
- symbol "Product safety class type BF";
- symbol "Danger of electric shock";
- IP classification;
- product weight;
- registration certificate number;
- barcode;
- Country of origin.


Fig. 11-3 - Marking of the plate of the main unit model Pro1000 on the back of the machine.

		Ultrasonic diagnostic device iLivTouch for non-	
		invasive determination of the degree of liver fibrosis	
Model Pro 3000 advanced search		SN _____	 _____
		I	
$\sim$ 100-240B, 50-60 Гц. 350 VA			 IPX0
			
	Wuxi Hisky Medical Technologies Co., Ltd Legal address: Room B401, 530 Plaza, University Science Park, Taihu International Science & Technology Park, 214135 Wuxi, Jiangsu Province, China Tel:+86-510-85387702. Web:www.hiskymedical.com		
Authorized representative of the manufacturer Limited Liability Company TESTSERVICE YUZHNY (LLC TESTSERVICE YUZHNY). Legal address: 143982, Russia, Moscow region, Balashikha, Yuzhnaya str. (Kuchino Mkr.), 9, floor 1, room, No 41. Tel.:7(499) 673-05-06 E-mail:regrzn_0506@mail.ru Registration certificate No. _____ dated _____			

Fig. 11-4 - Marking of the plate of the main unit model Pro3000 advanced search on the back of the machine.

11.7.3 Transport marking for delivery to the Russian Federation

The transport marking must comply with the requirements of GOST 14192 and must contain the following information (applied directly to the container and / or paper sticker):

- name and model of the medical device;
- full and abbreviated (if any) name of the manufacturer;
- manufacturer's postal address;
- full and abbreviated (if any) name of the authorized representative of the manufacturer;
- postal address of an authorized representative of the manufacturer;
- trademark of the manufacturer;

- "Serial number" symbol and related data;
- CE conformity mark;
- symbol date of manufacture and corresponding data (month, year);
- Overall dimensions of the package:
- number and date of registration certificate;
- storage conditions;
- transportation conditions;
- Symbol; Protect from moisture
- Symbol; Top
- Symbol; Fragile, handle carefully
- Symbol; Avoid sunlight
- Symbol; Stacking is prohibited;
- overall dimensions of the package with indication of the unit of measurement;
- Gross weight;
- Country of origin.

Ultrasonic diagnostic device iLivTouch for non-invasive determination of the degree of liver fibrosis

Model Pro 1000 **SN** _____  _____

 **HISKY MED**     

 **Wuxi Hisky Medical Technologies Co., Ltd** Legal address: Room B401, 530 Plaza, University Science Park, Taihu International Science & Technology Park, 214135 Wuxi, Jiangsu Province, China Tel:+86-510-85387702. Web:www.hiskymedical.com

Authorized representative of the manufacturer Limited Liability Company TESTSERVICE YUZHNY (LLC TESTSERVICE YUZHNY). Legal address: 143982, Russia, Moscow region, Balashikha, Yuzhnaya str. (Kuchino Mkr.), 9, floor 1. room. No 41. Tel.:7(499) 673-05-06 E-mail:regrzn_0506@mail.ru

Registration certificate No. _____ dated _____

Store at -5°C +40°C. relative humidity 0%-80%.

Transport at -20°C +50°C. relative humidity 30%-95%.

Gross weight: _____ kg. Overall dimension _____ x _____ x _____ cm

Made in China

Fig. 11-5 Layout of the transport labeling of the medical device model Pro 1000

**Ultrasonic diagnostic device iLivTouch for non-invasive
determination of the degree of liver fibrosis**

Model Pro 3000 advanced search SN

Wuxi Hisky Medical Technologies Co., Ltd Legal address: Room B401, 530
Plaza, University Science Park, Taihu International Science & Technology Park,
214135 Wuxi, Jiangsu Province, China Tel:+86-510-85387702.
Web:www.hiskymedical.com


Authorized representative of the manufacturer Limited Liability Company
TESTSERVICE YUZHNY (LLC TESTSERVICE YUZHNY). Legal address:
143982, Russia, Moscow region,
Balashikha, Yuzhnaya str. (Kuchino Mkr.), 9, floor 1, room, No 41.
Tel.:7(499) 673-05-06 E-mail:regrzn_0506@mail.ru
Registration certificate No. _____ dated _____
Store at -5°C +40°C, relative humidity 0%-80%.
Transport at -20°C +50°C, relative humidity 30%-95%.
Gross weight: _____ kg. Overall dimension ____x____x____cm
Made in China





Fig. 11-6. Layout of the transport labeling of the medical device model Pro 3000 advanced search

11.8.About symbols

The medical device uses the symbols contained in the technical documentation and the instruction manual, including Danger, Warning, Attention and Note. The meaning of these signs is described as follows:

Symbol	Value
 Danger	Indicates imminent hazards that could result in death or serious injury.
 Warning	Indicates a potential hazard that could result in death or serious injury.
 Attention	Indicates a potential hazard that could result in minor or moderate injury.

Symbol	Value
Note	Indicates instructions that operators must follow to avoid damage to the system and to ensure normal and efficient operation and function of the system.
	Global warning
	Working part type BF
OFF.	OFF for NETWORK management only
ON.	ON for NETWORK management only
	Equipotentiality
	Attention: Electric Shock Hazard
	Alternating current
IPX0	Not protected from water
IPX1	Protected against vertically falling water drops
IPX4	Splash-proof
IPX7	Protected against the effects of temporary immersion in water
IPX8	Protected against the effects of prolonged immersion in water
	Do not throw into trash. Use a separate collection for electrical and electronic equipment. Do not throw into trash. Use a separate
	See instruction manual/brochure NOTE on ME EQUIPMENT "Follow instructions for use"
	Manufacturer's declaration of conformity of the product with the applicable EEC directives, the European Notified Body is TÜV SÜD Product Service GmbH.
	Fuse
	Foot switch connector
	Top
	Fragile
	Protect from moisture
	Keep away from sunlight

Symbol	Value
	Stacking is prohibited;
	Ambient temperature for transportation: -20°C-50°C
	Ambient temperature for storage: -5°C-40°C
	Relative humidity during transportation: 30%-95%
	Relative humidity during storage: %
	Atmospheric pressure for transportation: 500hPa+1060hPa
	Atmospheric pressure for storage: 860hPa+1060hPa
	Main unit serial number
	Date of manufacture
	Manufacturer
	Authorized representative in European Community
	Additional information

11.9 Warranty and service life

HISKY MED guarantees products including the main unit and standard probes for eighteen (18) months from the date of shipment from the factory. We are responsible for troubleshooting and free replacement of damaged parts during the warranty period. Malfunctions specified in the warranty occur when the system is used in accordance with the Operation Manual.

The warranty period for products, including the main unit and standard probes, is eighteen (18) months from the date of shipment from the factory.

The expected life of the main unit and standard probes is 8 years.

The instruction manual must be enclosed in a polyethylene bag and placed in the shipping container together with the device.

They do not include the following:

1. Malfunctions, damage, defects caused by improper operation of the equipment;
2. Accidental or intentional damage;

3. Damage in transit;
4. Any product repaired or serviced by persons other than a distributor authorized by HISKY MED;
5. Man-made damage to the combination lock (we reserve the right to file a claim or legal action for damages caused by damage to the combination lock and related problems);
6. Defect due to force majeure such as flood, lightning, fire, etc.

Annex A: Acoustic power data

Sound power reporting table

Transducer model: XW-01

Working mode: Mode M

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.47	0.0003		0.0013		n/a
Index component value			0.0003	n/a	n/a	0.0013	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	0.75					
	P (MW)		0.024		0.024		n/a
	$P_{I \times I}$ (MW)		n/a		n/a		
	z_s (cm)		n/a				
	z_b (cm)					1.26	
	z_{MI} (cm)	1.26					
	$z_{P II, a}$ (cm)	1.26					
	f_{awf} (MHz)	2.57	2.57		2.57		n/a
Other information	p_{rr} (Hz)	39.10					
	s_{rr} (Hz)	n/a					
	n_{pps}	n/a					
	$I_{pa, a}$ at $z_{P II, a}$ (W/cm ²)	10.59					
	$I_{spta, a}$ at $z_{P II, a}$ or $z_{S II, a}$	0.24					
	I_{spta} at $z_{P II}$ or $z_{S II}$ (MW/cm ²)	0.30					
	p_r at $z_{S II}$ (MPa)	0.84					
Operational management conditions	Display depth (mm)	85	85	n/a	n/a	85	n/a
	Frequency (MHz)	2.5	2.5	n/a	n/a	2.5	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: B mode

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85		1.13		1.13	n/a
Index component value			1.13	1.13	n/a	1.13	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		109.9		109.9		n/a
	$P_{I \times I}$ (MW)		73.27		73.27		
	z_s (cm)			n/a			
	z_b (cm)					n/a	
	z_{MI} (cm)	4.38					
	$z_{PI,a}$ (cm)	4.38					
	f_{awf} (MHz)	3.25	3.25		3.25		n/a
Other information	p_{rr} (Hz)	3836.00					
	s_{rr} (Hz)	15.99					
	n_{pps}	1					
	$I_{pa,a}$ at $z_{PI,a}$ (W/cm ²)	72.26					
	$I_{spta,a}$ at $z_{PI,a}$ or $z_{SI,a}$	13.58					
	I_{spta} at z_{PI} or z_{SI} (MW/cm)	34.61					
	p_r at z_{SI} (MPa)	2.51					
Operational management conditions	Display focus (level)	9	9	9	n/a	9	n/a
	Display depth (mm)	120	120	120	n/a	120	n/a
	Operating frequency (MHz)	3.0	3.0	3.0	n/a	3.0	n/a
	Display focus number	1	1	1	n/a	1	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: Mode M

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85	0.040		0.11		n/a
Index component value			n/a	0.040	n/a	0.11	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		5.16		5.16		n/a
	$P_{1 \times 1}$ (MW)		n/a		n/a		
	z_s (cm)			3.02			
	z_b (cm)					4.18	
	z_{MI} (cm)	4.38					
	$z_{PII,a}$ (cm)	4.38					
	f_{awf} (MHz)	3.25	3.25		3.25		n/a
Other information	p_{rr} (Hz)	180.05					
	s_{rr} (Hz)	n/a					
	n_{pps}	n/a					
	$I_{pa,a}$ at $z_{PII,a}$ (W/cm ²)	72.26					
	$I_{spa,a}$ at $z_{PII,a}$ or $z_{SII,a}$	6.35					
	I_{spa} at z_{PII} or z_{SII} (MW/cm)	16.96					
	p_r at z_{SII} (MPa)	2.51					
Operational management conditions	Display focus (level)	9	n/a	9	n/a	9	n/a
	Display depth (mm)	120	n/a	120	n/a	120	n/a
	Operating frequency (MHz)	3.0	n/a	3.0	n/a	3.0	n/a
	Display focus number	1	n/a	1	n/a	1	n/a
	Scanning frequency (level)	High	n/a	High	n/a	High	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: Mode B+M

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85	1.10		1.17		n/a
Index component value			B:1,06 M: n/a	B:1,06 M: 0.040	n/a	B:1,06 M: 0.11	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		B:103,10 M:5,16		B:103,10 M:5,16		n/a
	$P_{1 \times 1}$ (MW)		B:68,73		B:68,73		
	z_s (cm)		3.02				
	z_b (cm)				4.18		
	z_{MI} (cm)	4.38					
	$z_{P_{II},a}$ (cm)	4.38					
	f_{awf} (MHz)	B:3,25	3.25		3.25		n/a
Other information	p_{rr} (Hz)	3600.96					
	s_{rr} (Hz)	15.00					
	n_{pps}	B:1					
	$I_{pa,a}$ at $z_{P_{II},a}$ (W/cm ²)	72.26					
	$I_{spta,a}$ at $z_{P_{II},a}$ or z_s (MW/cm ²)	19.09					
	I_{spta} at $z_{P_{II}}$ or $z_{S_{II}}$ (MW/cm ²)	49.43					
	p_r at $z_{S_{II}}$ (MPa)	2.51					
Operational management conditions	Display focus (level)	90	90	90	n/a	90	n/a
	Display depth (mm)	120	120	120	n/a	120	n/a
	Operating frequency (MHz)	3.0	3.0	3.0	n/a	3.0	n/a
	Display focus number	1	1	1	n/a	1	n/a
	Scanning frequency (level)	High	High	High	n/a	High	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: Mode B+CFM

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85	2.22	2.22	2.22	2.22	n/a
Index component value			B:1,60 C:0,62	B:1,60 C:0,62	n/a	B:1,60 C:0,62	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		B:155,00 C:68,48		B:155,00 C:68,48		n/a
	$P_{I \times I}$ (MW)		B:103,33 C:45,65		B:103,33 C:45,65		
	z_s (cm)			n/a			
	z_b (cm)					n/a	
	z_{MI} (cm)	4.38					
	$z_{PI,a}$ (cm)	4.38					
Other information	f_{awf} (MHz)	B:3,25	B:3,25 C:2,86		B:3,25 C:2,86		n/a
	prf (Hz)	3838.50					
	srf (Hz)	22.55					
	n_{pps}	B:1C:29					
	$I_{pa,a}$ at $z_{PI,a}$ (W/cm ²)	72.26					
	$I_{spta,a}$ at $z_{PI,a}$ or z_s	189.26					
	I_{spta} at z_{PI} or z_{SI} (MW/cm ²)	470.61					
Operational management conditions	p_r at z_{SI} (MPa)	2.51					
	Display focus (level)	9	9	9	n/a	9	n/a
	Display depth (mm)	120	120	120	n/a	120	n/a
	Operating frequency (MHz)	B:3,0 C:3,0	B:3,0 C:3,0	B:3,0 C:3,0	n/a	B:3,0 C:3,0	n/a
	Display focus number	1	1	1	n/a	1	n/a
	C PRF (kHz)	3.0	3.0	3.0	n/a	3.0	n/a
	Size (level)	15	15	15	n/a	15	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: PW mode

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.37		0.33		0.91	n/a
Index component value			n/a	0.33	n/a	0.91	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	0.63					
	P (MW)		44.20		44.20		n/a
	$P_{I \times I}$ (MW)		n/a		n/a		
	z_s (cm)			3.00			
	z_b (cm)					6.26	
	z_{MI} (cm)	6.50					
	$z_{P II, a}$ (cm)	6.50					
	f_{awf} (MHz)	2.88	2.88		2.88		n/a
Other information	p_{rr} (Hz)	5000.00					
	s_{rr} (Hz)	n/a					
	n_{pps}	n/a					
	$I_{pa, a}$ at $z_{P II, a}$ (W/cm ²)	18.38					
	$I_{spta, a}$ at $z_{P II, a}$ or z_s (MW/cm ²)	117.50					
	I_{spta} at $z_{P II}$ or $z_{S II}$ (MW/cm ²)	428.10					
	p_r at $z_{S II}$ (MPa)	1.20					
Operational management conditions	Display focus (level)	9	n/a	9	n/a	9	n/a
	Display depth (mm)	120	n/a	120	n/a	120	n/a
	Operating frequency (MHz)	3.0	n/a	3.0	n/a	3.0	n/a
	Display focus number	1	n/a	1	n/a	1	n/a
	PW PRF (kHz)	5.0	n/a	5.0	n/a	5.0	n/a
	Scanning frequency (level)	High	n/a	High	n/a	High	n/a
	Gateway (layer)	0	n/a	0	n/a	0	n/a
NOTE: N/A indicates that there is no corresponding intended use or no data is reported.							

Sound power reporting table

Transducer model: YD-CX

Working mode: Mode B+PW

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85	2.59		3.17		n/a
Index component value			B:2,26 PW: n/a	B:2,26 PW:0,33	n/a	B:2,26 PW:0,91	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		B:218,65 PW:44,20		B:218,65 PW:44,20		n/a
	$P_{1 \times 1}$ (MW)		B:145,77		B:145,77		
	z_s (cm)			3.00			
	z_b (cm)					6.26	
	z_{MI} (cm)	4.38					
	$z_{PII,a}$ (cm)	4.38					
	f_{awf} (MHz)	B:3,25	B:3,25 PW:2,88		B:3,25 PW:2,88		n/a
Other information	p_{rr} (Hz)	3618.20					
	s_{rr} (Hz)	31.81					
	n_{pps}	1					
	$I_{pa,a}$ at $z_{PII,a}$ (W/cm ²)	72.26					
	$I_{spa,a}$ at $z_{PII,a}$ or $z_{SII,a}$	144.52					
	I_{spa} at z_{PII} or z_{SII} (MW/cm)	496.96					
	p_r at z_{SII} (MPa)	2.51					
Operational management conditions	Display focus (level)	9	9	9	n/a	9	n/a
	Display depth (mm)	120	120	120	n/a	120	n/a
	Operating frequency (MHz)	B:3,0 PW:3,0	B:3,0 PW:3,0	B:3,0 PW:3,0	n/a	B:3,0 PW:3,0	n/a
	Display focus number	1	1	1	n/a	1	n/a
	PW PRF (kHz)	5.0	5.0	5.0	n/a	5.0	n/a
	Scanning frequency (level)	High	High	High	n/a	High	n/a
	Gateway (layer)	0	0	0	n/a	0	n/a

NOTE: N/A indicates that there is no corresponding intended use or no data is reported.

Sound power reporting table

Transducer model: YD-C X

Working mode: Mode B+CFM+PW

Index label		MI	TIS		TIB		TIC
			On a surface	Under surface	On a surface	Under surface	
Maximum index value		0.85		1.11		1.69	n/a
Index component value			B:0,56 C:0,22 PW: n/a	B:0,56 C:0,22 PW: 0.33	n/a	B:0,56 C:0,22 PW: 0.91	
Acoustic parameters	$p_{r,a}$ at Z_{MI} (MPa)	1.53					
	P (MW)		B:53,96 C:23,84 PW:44,20		B:53,96 C:23,84 PW:44,20		n/a
	$P_{I \times I}$ (MW)		B:35,97 C:15,89		B:35,97 C:15,89		
	z_s (cm)			3.00			
	z_b (cm)					6.26	
	z_{MI} (cm)	4.38					
	$z_{P II, a}$ (cm)	4.38					
Other information	f_{awf} (MHz)	B:3,25	B:3,25 C:2,86 PW:2,88		B:3,25 C:2,86 PW:2,88		n/a
	prf (Hz)	3618.20					
	srr (Hz)	7.85					
	n_{pps}	1					
	$I_{pa, a}$ at $z_{P II, a}$ (W/cm ²)	72.26					
	$I_{spa, a}$ at $z_{P II, a}$ or $z_{S II, a}$	183.39					
	I_{spa} at $z_{P II}$ or $z_{S II}$ (MW/cm ²)	591.92					
Operational management conditions	p_r at $z_{S II}$ (MPa)	2.51					
	Display focus (level)	9	9	9	n/a	9	n/a
	Display depth (mm)	120	120	120	n/a	120	n/a
	Operating frequency (MHz)	B:3,0 C:3,0 PW:3,0	B:3,0 C:3,0 PW:3,0	B:3,0 C:3,0 PW:3,0	n/a	B:3,0 C:3,0 PW:3,0	n/a
	Display focus number	1	1	1	n/a	1	n/a
	CFM PRF (kHz)	3.0	3.0	3.0	n/a	3.0	n/a
	Size (level)	15	15	15	n/a	15	n/a
	PW PRF (kHz)	5.0	5.0	5.0	n/a	5.0	n/a
	Scanning frequency (level)	High	High	High	n/a	High	n/a
	Gateway (layer)	0	0	0	n/a	0	n/a

NOTE: N/A indicates that there is no corresponding intended use or no data is reported.

Appendix B: EMC Declaration of Conformity

Warning:

1. Use of this equipment near or stacked with other equipment should be avoided as it may cause improper operation. If such use is necessary, this system and other equipment should be observed to ensure that they are operating correctly.
2. Use of accessories, probes, and cables other than those specified or supplied by the manufacturer of this equipment may result in increased electromagnetic emissions or reduced immunity of this equipment, and may result in improper operation.
3. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used at least 30 cm (12 inches) away from any part of the system, including cables specified by the manufacturer. Otherwise, it may degrade the performance of this equipment.

Note:

1. The EMISSION characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If the equipment is used in a residential environment (which typically requires CISPR 11 Class B), it may not provide adequate protection for RF communications services. The user may be required to take risk mitigation measures such as relocating or reorienting the equipment.
2. Use cables and other accessories that meet the requirements from the following list:

Port No.	Name:	Type*	Maximum cable length >3 m	Shielded cable
1	Net	Alternating current	<3m	Unshielded
2	USB-port	Input/Output	0 m	—
3	Fibrosis probe socket	Input/Output	<3m	Shielded
4	Imaging probe socket	Input/Output	<3m	Unshielded
5	Foot switch port	Input/Output	<3m	Unshielded

This product complies with the IEC 60601-1-2 EMC standard: 2014. Compliance with

each EMISSION and IMMUNITY standard or test specified in this supplementary standard is shown in the EMI (electromagnetic emission) and EMI (electromagnetic susceptibility) tables below.

Correspondence table for EMI and EMW

Table 1 - Emission

Phenomenon	Correspondence	Electromagnetic environment
RF radiation	CISPR 11 Group 1, Class A	Professional healthcare environment
Harmonic distortion	IEC 61000-3-2 Class A	Professional healthcare environment
Voltage fluctuations and flicker	IEC 61000-3-3 Correspondence	Professional healthcare environment

EMC Compliance Table

Table 2 - Body port

Phenomenon	Basic EMC standard	Susceptibility test levels
		Professional healthcare environment
Electrostatic discharge	IEC 61000-4-2	contact ± 8 kV ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air
Radiated RF electromagnetic field	IEC 61000-4-3	3 V/m 80MHz-2.7GHz 80% AM, 2Hz or 1kHz
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	See table 3
Magnetic fields of rated power frequency	IEC 61000-4-8	30A/m 50Hz or 60Hz

Table 3 - Proximity fields from RF wireless communications equipment

Test frequency (MHz)	Range (MHz)	Susceptibility test levels
		Professional healthcare environment
385	380-390	Pulse modulation 18 Hz, 27 V/m
450	430-470	Frequency modulation ± 5 kHz sine deviation 1 kHz, 28 V/m
710	704-787	Pulse modulation 217 Hz, 9 V/m
745		
780		
810	800-960	Pulse modulation 18 Hz, 28 V/m
870		
930		
1720	1700-1990	Pulse modulation 217 Hz, 28 V/m
1845		
1970		
2450	2400-2570	Pulse modulation 217 Hz, 28 V/m
5240	5100-5800	Pulse modulation 217 Hz, 9 V/m
5500		
5785		

Table 4 - AC Power Input Port

Phenomenon	Basic EMC standard	Susceptibility test levels
		Professional healthcare environment
Fast electrical transients/bursts	IEC 61000-4-4	± 0.5 kV, ± 1 kV, ± 2 kV Repetition frequency 100 kHz
Bursts Line voltage	IEC 61000-4-5	± 0.5 kV, ± 1 kV
Bursts Phase voltage	IEC 61000-4-5	± 0.5 kV, ± 1 kV, ± 2 kV
Conducted interference caused by RF fields	IEC 61000-4-6	3V, 0.15MHz-80MHz 6 V on ISM bands from 0.15 MHz to 80 MHz 80%AM at 1 kHz and 2 Hz
Voltage drops	IEC 61000-4-11	0% U_T ; 0.5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270° и 315°
		0% U_T ; 1 cycle and 70% U_T ; 25/30 cycle Single phases at 0°
Voltage interruptions	IEC 61000-4-11	0% U_T ; 250/300 cycle At 0°, 180°

Application C: List of international and national standards

International Standards

1. IEC60601-1:2005/A1:2012 Medical Electrical Products - Part 1: General requirements for basic safety and basic functional characteristics";
2. IEC60601-1-2:2014 Medical Electrical Products - Part 1-2: General requirements for basic safety and basic performance - Parallel standard: Electromagnetic compatibility - Requirements and tests;
3. IEC60601-2-37:2007/A1:2015 Medical electrical products - Part 2-37: Additional requirements for basic safety and basic functional characteristics of ultrasonic medical equipment for diagnostics and control;
4. IEC 60601-1-6:2010/A1:2013 Medical electrical products - Part 1-6: General requirements for basic safety and basic performance - Parallel standard: Practicality;
5. IEC 62304:2006/A1:2015 Medical device software - Software life cycle processes;
6. IEC 62366-1:2015/A1:2016 Medical devices - Part 1: Applying usability design to medical devices;
7. ISO 10993-1:2018 Biological evaluation of medical devices -- Part 1: Assessment and testing as part of the risk management process;
8. ISO 10993-5:2009 Biological evaluation of medical devices -- Part 5: Tests for cytotoxicity in laboratory conditions (ISO 10993-5:2009);

9. ISO 10993-10:2010 Biological evaluation of medical devices -- Part 10: Tests for skin sensitization;
10. ISO 14971:2019 Medical Devices - Applying risk management to medical devices;
11. ISO 13485:2016 Medical devices - Quality management systems - Regulatory requirements;
12. ISO 15223-1:2016 Medical products. Symbols used for marking on medical devices, labels and accompanying documentation - General requirements;
13. EN 1041:2008/A1:2013 Information provided by the medical device manufacturer.

National standards (compliance with Russian standards)

1. GOST R 50444-2020 "Medical devices, apparatus and equipment. General technical conditions";
2. GOST R MEC 60601-1-2010 "Medical electrical device. Part 1. General safety requirements, taking into account the basic functional characteristics.
3. GOST R MEC 60601-1-2014 "Medical electrical device. Part 1-6. General safety requirements, taking into account the basic functional characteristics";
4. GOST R IEC 62304-2013 "Medical products. Software Life cycle processes";
5. GOST R IEC 62366-2013 "Design of medical devices taking into account serviceability";
6. GOST ISO 14971-2011 "Medical devices - Application of risk management to medical devices";
7. GOST R ISO 15223-1-2020 "Medical devices. Symbols used for marking on medical devices, labels and accompanying documentation."

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