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EmbryoSlide+™ culture dish

40 EmbryoSlide+ culture dishes with lids
1 dish per bag

Indications for use: The EmbryoSlide+ culture dish is intended for preparing, storing and transferring human embryos. The EmbryoSlide+ culture dish must be used together with the EmbryoScope+ incubator.

1-cell MEA ≥ 80% expanded blastocysts at 96 hrs.

STERILE Sterilized using irradiation with SAL 10⁶ if packaging sealed and undamaged.

REF 16450

Commercial Name

LOT P-P8-2110-1

Batch Code

2021-10-01

Manufactured Date

2025-10-01

Expiry

UDI

(01)5712714664509(10)PP821101(17)261001

Do not use if package is damaged

Do not reuse

Consult instructions for use

Store at room temperature

Vitrolife A/S
Jens Juuls Vej 20
8260 Viby J, Denmark

2460

Caution: Federal law restricts this device to sale by or on the order of a physician or a practitioner trained in its use.

Vitrolife

For contact

Americas

Tel: + 1-866-848 7687

support.embryoscope.us@vitrolife.com

Europe, Middle East, Asia, Pacific, Africa

Tel: +45 8944 9500

support.embryoscope@vitrolife.com

Vitrolife A/S, Jens Juuls Vej 20, 8260 Viby J, Denmark

Vitrolife Sweden AB, Box 9080 (Gustaf Werners gata 2), SE-400 92 Göteborg, Sweden

Vitrolife Inc., 3601 South Inca Street, Englewood, CO 80110, USA

Vitrolife Pty. Ltd., 55 Flemington Rd, North Melbourne VIC 3051, Australia

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EN

Intended use

The EmbryoSlide+ culture dish is an accessory for embryo storage related to the EmbryoScope providing an environment with controlled temperature and CO₂ for the development of embryos at or near body temperature.

Indications for use

The EmbryoSlide+ culture dish is intended for preparing, storing and handling human embryos. The EmbryoSlide+ culture dish must be used together with the EmbryoScope incubator.

Restrictions on use

EmbryoSlide+ is only intended for use in the EmbryoScope+ incubator and must be handled by trained personnel according to the instructions in this insert.

This EmbryoSlide+ culture dishes and lids are sterile. The pouches must only be opened in a sterile laminar flow hood and the dish must be covered by the lid when in use.

Precautions and warnings

WARNING: EmbryoSlide+ is intended for single use only and MAY NOT be re-used. Any attempt by the user to clean and re-use the dish may result in contamination with microorganisms or other risks of device failure.

WARNING: To avoid contamination with microorganisms, always place the culture dish in a sterile laminar flow hood while handling and generally handling the dish.

WARNING: Always label the culture dish appropriately and indicate the label when the embryo is transferred either to a new device or to the patient.

CAUTION: Follow aseptic technique to avoid the risk of infection of a physician or a practitioner involved in its use.

CAUTION: Do not use the mixing wells for incubating embryos as no oxygen is supplied from these wells.

CAUTION: The reservoir and mixing wells MUST always be covered by a confluent oil layer of 1.5 ml, IVF-grade oil. This prevents the medium from evaporating during incubation.

CAUTION: Evaporation of the medium can change osmolality, which may affect embryo development.

CAUTION: Always place the dish in a sterile laminar flow hood while loading it to keep both hands free and reduce the risk of spilling any oil or medium.

CAUTION: If ANY oil or medium is accidentally spilled from the dish while loading or handling it, the embryo MUST be transferred to a new dish to prevent them from being adversely affected by the spillage.

CAUTION: In case any oil is spilled on the culture dish outside of the reservoir, gently wipe off the oil by using a lint-free tissue paper. If any oil is present on the lid, the lid must be replaced.

CAUTION: If any bubbles are present in the well or oil layer after preparation, they may block the laminar flow and compromise image clarity. Oil-free bubbles may damage the embryo from the bottom of the reservoir.

CAUTION: Any bubbles must be removed carefully and reactively quickly to avoid evaporation of medium.

CAUTION: Take care not to remove any medium when removing the bubbles.

Preparing the culture dish

Prepare one culture dish at a time to minimize handling time. Work on a non-heated workbench and use cold medium and oil.

1. Remove culture dish from pouch in a sterile laminar flow hood.
2. Fill an incubator with medium using a micropipette (max. diameter 200 µm). One filling of the tip will suffice to completely fill eight microwells. To avoid bubbles, slightly overfill the microwell to avoid creating a central depression in the surface of the medium. The tip of the pipette must touch the side of the microwell during filling.
3. Use a standard pipette to fill an additional 180 µl of medium into the wells and the two reservoirs. Rotate the tip over the wells while releasing medium. Make sure to fill the reservoirs completely, including the spilling zone.
4. Fill each mixing well with max. 30 µl, and min. 25 µl, of medium.
5. Load 1.5 ml of oil verified for IVF into the reservoir. Make sure that all culture wells and mixing wells are covered with oil.
6. Push up larger bubbles with a micropipette and then remove them.
7. Cover with lid and equilibrate overnight.
8. Identify and remove any bubbles under a stereo microscope.
9. Place the dish in the incubator and label the lid on the dish. The bubble must be covered so that it is possible to read the label when the dish is correctly inserted in the incubator. Make sure that the label is smooth and not as best as possible to avoid wrinkles in the paper. If the bubble is in any way damaged or not smoothly applied, the label cannot be read.
10. Load embryos into the centre of the microwells.
11. Place the dish in the EmbryoScope incubator.

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