Medtronic

Further. Together

Medtronic Inc.

710 Medtronic Parkway


Minneapolis, MN 55432

USA

Ashley Heerema
 Ashley Heerema
 Sr. Regulatory Affairs Specialist
 Date *Dec 10, 2019*
 Medtronic, Inc.

Addendum to Device Manual

1. Name of medical product	Pacemaker implantable single-chamber Azure XT SR MRI SureScan, model W2SR01
2.1. Prescription of a medical product	Medical device is a multiprogrammable cardiac device, which monitors and regulates the patient's heart rate, by performing a one-chamber rate-responsive adaptive stimulation with bradycardia.
2.2 Description of the principles underlying the work of the medical products and their features	Feature summary Pacing features Auto-adjusting sensitivity – This feature automatically adjusts the sensitivity thresholds after specific paced events and sensed events occur. Automatic polarity configuration – This device uses lead impedance measurements to automatically configure pacing and sensing polarities during Implant Detection. Conducted AF Response – This feature regularizes the ventricular rhythm during conducted AT/AF by modifying the pacing rate on a beat-by-beat basis to closely match the patient's average ventricular response. MRI SureScan – This feature allows patients with an implanted MRI SureScan system, including the device and lead, to have a safe MRI procedure if the requirements provided in the MRI technical manual are followed.

Cindi A Bigelbach
12/10/2019

Rate Profile Optimization – The goal of Rate Profile Optimization is to ensure that the rate response remains appropriate for the full range of patient activities. This feature monitors the patient's daily and monthly sensor rate profiles and adjusts the rate response curves over time to achieve a prescribed target rate profile.

Rate-responsive pacing – This feature varies the pacing rate in response to the patient's physical motion as detected by the activity sensor of the device.

RV Capture Management – This feature monitors the right ventricular pacing threshold with daily pacing threshold searches and, if programmed to do so, adjusts the RV pacing amplitude toward a target amplitude.

Sleep feature – This feature causes the device to pace at a slower rate during a programmed sleep period.

Ventricular Rate Stabilization (VRS) – This feature adjusts the pacing rate dynamically to eliminate the long pause that typically follows a premature ventricular contraction (PVC).

Monitoring and follow-up features

Cardiac Compass Trends – This feature provides a Cardiac Compass Trends report that shows an overview of the patient's condition, with graphs that display long-term trends in heart rhythm over the last 14 months. The report also includes the OptiVol 2.0 fluid trend data.

Medtronic CareAlert Monitoring – If the device identifies any programmed or automatic CareAlert conditions, this feature sends a wireless alert signal to the patient monitor (if available). The patient monitor then transmits the CareAlert Event data to the Medtronic CareLink Network. If configured to do so, the Medtronic CareLink Network then sends an alert notification to the clinic.

Episode data and EGM storage – The system provides an arrhythmia episode log that enables you to view the summary and detailed diagnostic data quickly, including stored EGM, for the selected arrhythmia episode.

Flashback memory – This diagnostic feature records intervals that occur immediately prior to tachyarrhythmia episodes or the most recent interrogation and plots the interval data over time.

Heart Failure Management Report – This report provides an overview of the patient's condition over the short and long term, with a focus on heart failure management. The report includes graphs that show OptiVol 2.0 fluid trends and trends related to heart failure over the last 14 months.

Holter telemetry – This function allows the implanted device to transmit an EGM with marker telemetry continuously for up to 46 hours, regardless of the use of the programming head.

Implant Detection – Implant Detection is a 30 min period, beginning when the device is placed in the surgical pocket. During this period, the device verifies lead connection by measuring lead impedance. When the Implant Detection period is completed, various automatic features and diagnostics are activated.

Lead Monitor – This feature measures lead impedances during the life of the implanted device and controls automatic configuration of lead polarities at implant. If Lead Monitor is programmed to Adaptive, the device automatically switches bipolar pacing and sensing to unipolar pacing and sensing if the integrity of a bipolar lead is compromised.

OptiVol 2.0 fluid trends – This feature provides the capability to monitor the following trends:

- The Thoracic Impedance trend plots thoracic impedance for up to 14 months.
- The OptiVol 2.0 Fluid Index trend plots the accumulated differences between the Daily Impedance and Reference Impedance values. Possible fluid accumulation in the patient's thoracic cavity exists when the OptiVol 2.0 Fluid Index exceeds the OptiVol Threshold.

Rate Histograms report – This report shows heart rate range distributions for the patient.

Data security

Medtronic has designed safeguards to protect patient information and device data for the Astra XT SR MRI SureScan Model X2SR01 device.

Bluetooth communication system – The device shows its availability through Bluetooth communication. Critical data accepted or sent through the Bluetooth communication from the device is encrypted by the device before it is sent over the Bluetooth channel. The device responds only to authorized commands.

Inductive telemetry communication system – The Medtronic inductive telemetry communication system is used with the clinician programmer to interrogate and program the device. It can also be used to interrogate the device for remote monitoring, if available. This system uses short-range communication that protects patient information and device data.

Pacing mode information

Pacemaker modes are described using the NBG code. The five-letter NBG code, named after The North American Society of Pacing and Electrophysiology (NASPE) and the British Pacing and Electrophysiology Group (BPEG), describes the operation of implantable pulse generators. The NBG code, which supersedes the ICHD Code, is described in Table 1.

	Position:	I	II	III	IV	V
	Category:	Chamber(s) Paced	Chamber(s) Sensed	Response to Sensing	Rate Modula- tion	Multisite Pac- ing ^a
		O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None T = Triggered I = Inhibited D = Dual (T + I)	O = None R = Rate mod- ulation	O = None A = Atrium V = Ventricle D = Dual (A + V)
	Manufacturers' designation only:	S = Single ^b (A or V)	S = Single ^b (A or V)			
^a Medtronic devices do not use the Multisite Pacing code.						
^b The programmer displays A or V (not S) for chambers paced and sensed.						
3. Conditions for use	In medical and preventive treatment facility. The pacemaker should be used only by or under the supervision of physicians trained.					
4. Composition and list of parts and accessories*	Implantable pacemaker single chamber Azure XT SR MRI SureScan, model W2SR01 consisting of: - implantable pacemaker single chamber Azure XT SR MRI SureScan, model W2SR01; - torque wrench -Product documentation includes - Device manual; Limited Warranty; Medical procedures and EMI Warnings and Precautions manual; MRI Technical manual ; Explanation of the symbols IPG and CRT-P devices; Stickers; Registration card; Medtronic Web Card; CD-ROM. * Appendix to the Device Manual isn't included in package content, provided to customer separately on the territory of Russian Federation, SW030 Programmer software guide and Reference Manual isn't included in package content, visit http://manuals.medtronic.com/					
5. Technical details of medical product	Pacemaker implantable single-chamber Azure XT SR MRI SureScan, model W2SR01:					
	Volume ^a			12.25 cm ³		
	Mass			22.5 g		
	H x W x D ^b			42.6 mm x 50.8 mm x 7.4 mm		
	Radiopaque ID ^c			RNA		

Medtronic identifier	
Surface area of titanium device can	33.48 cm ²
Materials in contact with human tissue ^d	Titanium, polyurethane, silicone rubber
Battery	Lithium-hybrid CFx silver vanadium oxide

^a Volume with connector holes unplugged.





^b Grommets may protrude slightly beyond the can surface.

^c The radiopaque ID, which includes a Medtronic-identifier symbol, can be viewed in a fluoroscopic image of the device.

^d These materials have been successfully tested for the ability to avoid biological incompatibility. The device does not produce an injurious temperature in the surrounding tissue during normal operation.

	Height (mm)	Width (mm)	Case thickness (mm)	Connector Thickness (mm)	Mass (g)	Volume (cm ³)
Single Chamber Mode	42.6	50.8	7.4	9.9	22.5	12.25
Tolerance	+/-10%	+/-10%	+/-10%	+/-10%	+/-10%	+/-10%


Radiopaq ID drawing and location:

SECTION E-E

All measurements and tolerances for radiopaque drawing is in inches.
Unless otherwise specified, default tolerance is ± 0.006 inches


Torque wrench (included in the package with ICD) – Use only the torque wrench supplied with the device. The torque wrench is designed to prevent damage to the device from overtightening a setscrew. Other torque wrenches (for example, a blue-handled or right-angled hex wrench) have torque capabilities greater than the lead connector can tolerate.

*

	The torque is limited to 12 to 17 in-oz. during first clockwise revolution; it will click when torque value in this range is achieved Mass: 8 g $\pm$ 0,5 g *Tolerances of dimensions are $\pm$ 5%, unless otherwise specified Software Model of SW - SW030 Version – 8.0 and upper Safety Class - C
6. The possibility and ways of integrating with other medical products	The main components of the system are the IPG, SureScan lead(s), and physician programmer. The patient monitor (available with Astra and Azure models) or wireless monitoring systems (available with Azure models only) are optional system components allowing patients to communicate IPG data through the CareLink Network to clinicians for review (Figure 1). The Patient Assistant is also an optional part used for marking symptoms. This clinical evaluation report (CER) covers the IPG component and the programmer application software, the remaining system components are covered by separate CERs. Astra and Azure devices sense electrical activity of the patient's heart using electrodes in the implanted leads and then analyze heart rhythms based on selectable therapy parameters and provide therapy as appropriate. These devices also provide diagnostic and monitoring information for physician-directed patient care and system evaluation. These devices are Magnetic Resonance (MR) conditional and allow Magnetic Resonance Imaging (MRI) under certain conditions.

Figure 1: Astra and Azure System Architecture

Only the physician programmer, using inductive telemetry communication system is capable of programming the Astra and Azure IPGs. Patients may use the optional Medtronic patient monitor to gather information from their implanted devices and to send information to their physicians through the Medtronic CareLink Network (patient-initiated transmissions). The alternative, optional, wireless Bluetooth Low Energy (BTLE) mobile application (not subject to medical device regulation) and BTLE enabled Medtronic patient monitor (available for Azure models only) allow device-initiated (CareAlerts) and remote scheduled transmissions without patient initiation. The IPG will only respond to a restricted set of commands, such as interrogation commands. These commands do not allow therapy parameters to be modified. Data sent from the implanted device through the wireless communication system will be encrypted in the IPG and will pass through the mobile app to the CareLink server where it is decrypted. The transmissions between the IPG and CareLink server are encrypted end-to-end.

Description of Accessories and product connected to other devices

The Astra and Azure MRI Systems may be comprised of the following components:

- CareLink Programmer Model 2090
- CareLink Encore Programmer Model 29901
- Conexus Activator Model 27901
- MyCareLink Monitor Model 24952
- MyCarelink Smart, Model 25000 (Astra only)
- 2020B CareLink Express Monitor
- MyCareLink Heart Model 31302 and Model 31301 (Azure only)
- Commercially available MRI SureScan Right Atrial Pacing Leads¹
- Commercially available MRI SureScan Left Ventricular Pacing Leads²
- Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected)
- The same commercially available implant support instruments and accessories used with the predecessor Advisa DR device systems.

Commercially Available Accessories

The Azure/Astra system will be compatible with the following commercially available system components:

- Model 2290 pacing system analyzer
- Model 9466 Patient Magnet
- Model DR220 Holter Monitor
- Model PA96000 Patient Assistant

¹ The connector ports are designed according to the IS-1 standard ISO 5841-3 or the IS4 / DF4 standard ISO 27186. These ports are compatible with the lead configurations as specified in the product labeling.


	Lead Model	Brand Name	
	5086MRI	CapSureFix MRI SureScan Pace/Sense Lead	
	5076	CapSureFix Novus MRI SureScan Pace/Sense Lead	
	4076	CapSureFix Novus MRI SureScan Lead	
	4074/4574	CapSure Sense MRI SureScan Lead	
	3830	SelectSecure Model Pacing Lead	
	5554	CapSure Z Novus Model 5554 pacing lead	
	5054	CapSure Z Novus Model 5054 pacing lead	


7. EMC Information	Electromagnetic emission		
	Medical device are intended for use in the electromagnetic environment defined below. Buyer or user medical device should ensure its use in the specified environment.		
	Electromagnetic emission test	Conformity	Electromagnetic environment - instructions
	Radio interference according to CISPR 11	Group 1	Medical device should emit electromagnetic energy to perform the main function.
	Radio interference to CISPR 11	Class B	Medical device are suitable for use in all locations.
	Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use	

		Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use			
		Noise immunity				
		Immunity test	The test level of IEC 60601	Level of compliance	Electromagnetic environment - instructions	
		Electrostatic Discharge (ESD) in accordance with IEC 61000-4-2	±6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%	
		Magnetic field of industrial frequency according to IEC 61000-4-8	3 A/m	0,3 A/m	The magnetic fields of the industrial frequency must be measured at the designated installation site to ensure that the field strength is sufficiently low	
		Immunity test	The test level of IEC 60601	The actual level of noise immunity	Level of compliance	
		Conductive noise induced by radio-frequency electromagnetic fields in	3 V (rms) in the band from 150 kHz to 80 MHz outside the frequencies allocated for the	1.7 V (rms value)	1 V (rms value)	


		accordance with IEC 61000-4-6	LPM of HF devices		
		The emitted radio-frequency electromagnetic field according to IEC 61000-4-3	10 V (RMS) in the band from 150 kHz to 80 MHz at the frequencies allocated for the LPM of high-frequency devices	2.3 V (rms value)	1 V (rms value)
				13 V/m	10 V/m
		Noise immunity			
		The test for Noise immunity	Test Level according to IEC 60601	Level Conformity	Electromagnetic environment - instructions
		Electrostatic Discharge (ESR) according to IEC 61000-4-2	± 6 kV – contact discharge ± 8 kV – air discharge	± 6 kV – contact discharge ± 8 kV – air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%.
		A magnetic field Industrial Frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of the industrial frequency must correspond to the typical conditions of a commercial or hospital environment.
		Radio frequency Electromagnetic field IEC 61000-4-3	3 V / m In the band from 80 MHz to 2.5 GHz	3 V/m	However, the readings of the weights built into the bed do not have a critical value for the operation of a medical device: D = $1,2\sqrt{P}$ in the band from 26 to 800 MHz D = $2,3\sqrt{P}$ in the band from 800 MHz to 2.5 GHz,

					Where P is the nominal maximum output transmitter power in Watts (W) set by the manufacturer, and d - Recommended spatial clearance in meters (m). The field strength in the propagation of radio waves from stationary radio transmitters from the results of observations of the electromagnetic situation should be lower than the level of compliance in each frequency band. Interference may occur near equipment marked with a sign.  This device may not cause harmful interference, and this to accept any interference received, including interference that may cause undesired operation of the device. This device must not cause harmful interference to radio stations operating in the 400.150-406.000 MHz band and used for meteorological purposes, for meteorological satellites and systems of Earth exploration satellites. This device must accept any interference received, including interference that may cause undesired operation.	
--	--	--	--	--	---	--

	NOTE 1 Higher frequency bands are applied at 80 MHz and 800 MHz. NOTE 2: This manual is not applicable in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.
8. Materials used for production of medical device	Implantable pacemaker Titanium grade 1, Silicone rubber, Polyurethane Torque wrench Polyetherimide, Corrosion resistant steel
9. Design and appearance, photo or image	Pacemaker implantable single-chamber Azure XT SR MRI SureScan, model W2SR01 consisting of: Pacemaker implantable single-chamber Azure XT SR MRI SureScan, model W2SR01  Torque wrench


Product documentation


10. Information on labelling, packaging

Sterile barrier systems constructed of uncoated PETG (Glidex) trays, sealed to Oliver 18B coated Tyvek® 1073B porous lid stock and sterilized with ethylene oxide gas. This sterile barrier material configuration is common to many CRDM device models across different product families. Much of the testing that supports the use of packaging materials and sterile package assemblies is specific to the packaging materials and assemblies, not specific products. Product specific evidence for sterilization qualification, package performance qualification, labeling, and product specific requirements can be found in the DHF for each individual product.


CRDM Sterile Barrier systems constructed of uncoated PETG trays, sealed to Oliver 18B coated Tyvek® 1073B porous lid stock and sterilized with ethylene oxide gas include the IPG/ICD/CRT Device Package and the Dual Entry Leads Package. Images of these packages follow:


Package dimensions:
1. Sterile package lid (in inches)
Inner package lid:


Outer package lid:

Sterile Package dimensions:

Outer package Drawing with Dimensions:

***Tolerances of dimensions for each package are $\pm 5\%$, unless otherwise specified**

Labeling of medical devices packages:

Pacemaker implantable single-chamber Azure XT SR MRI SureScan, model W2SR01 (shelf box and single package labels)

Medtronic **Bluetooth**

Azure™ XT SR
MRI SureScan™

W2SR01
SN WWW000000W
2000-00-00
2000-00-00
REF W2SR01

VVIR IS-1 V SureScan MR

1x

VVI 60 min⁻¹
0.9 mV
3.5 V
0.4 ms
AUTO

Adaptive

STERILE EO -18 °C 0 °F 55 °C 131 °F www.medtronic.com/manuals

FCC ID: LFSBLEIMPLANT
IC: 3408D-BLEIMPLANT

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands
Manufactured in: Juncos, Puerto Rico, USA
Product of USA
www.medtronic.com

CE0123/1856

Medtronic **Bluetooth**

Azure™ XT SR MRI
SureScan™

W2SR01
SN WWW000000W
2000-00-00
2000-00-00
REF W2SR01

VVIR IS-1 V SureScan MR

VVI 60 min⁻¹
0.9 mV
3.5 V
0.4 ms
AUTO

Adaptive

STERILE EO -18 °C 0 °F 55 °C 131 °F www.medtronic.com/manuals

FCC ID: LFSBLEIMPLANT
IC: 3408D-BLEIMPLANT

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands
Manufactured in: Juncos, Puerto Rico, USA
www.medtronic.com

CE0123/1856

M970311A001 B - M967972A001 A

2000-00-00
YYYY-MM-DD

Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W
Medtronic W2SR01	 	Azure™ XT SR MRI SureScan™ SN WWW000000W

Medtronic

Azure™ XT SR MRI SureScan™ W2SR01

Implantation
date:

2000-00-00
YYYY-MM-DD

REF W2SR01

SN WWW000000W

Medtronic

Azure™ XT SR MRI SureScan™ W2SR01

Implantation
date:

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

SN WWW000000W

Medtronic

Azure™ XT SR MRI SureScan™ W2SR01

Implantation
date:

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

SN WWW000000W

Medtronic

VVIR

1x IS-1

Azure™ XT SR MRI SureScan™ W2SR01

• BlueSync™ cardiac monitoring system

REF

W2SR01

SN

WWW000000W

2000-00-00
YYYY-MM-DD

PIN: W000WWW00W

M970311A001 B - M967972A001 A

Medtronic

VVIR

1x IS-1

Azure™ XT SR MRI SureScan™ W2SR01

(01)00000000000000(17)000000(21)WWW000000W

REF

W2SR01

SN

WWW000000W

2000-00-00
YYYY-MM-DD

M970311A001 B - M967972A001 A

Sticker for final package will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation):

	Product name (in Russian)	
	Manufacturer (name)	UPN (unique product number)
	Country of origin	CFN n (customer facing number)

See IFU

RC # _____ date of issue _____

The below list of symbols applies to various products. Refer to the product package labels to see which symbols apply to this product.

Products and Systems incorporating radio frequency transmitters or intentionally using radio frequency electromagnetic energy for diagnosis or treatment shall be marked in the form of a sign of non-ionizing radiation

This conformity mark is placed on the products for which the mandatory certification procedure was carried out

Packing materials:

Package part	Material
Trays	PETG
Lids	Tyvek® 1073B
Label on primary package	Inkjet paper
Print ink blue	Blue ink
Print ink black	Black ink

11. Methods and conditions of sterilization and time

Sterilized by EO. The Azure and Astra Implantable Pulse Generator (IPG) product family is qualified for 1X – 5X sterilization in the 30-minute 100% EO sterilization process.

maintaining sterility (for sterile medical devices produced)	
12. Requirements for the protection of the environment in the application of medical products	Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See http://recycling.medtronic.com for instructions on proper disposal of this product. Consider the following information related to device explant and disposal: ● To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. ● Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. ● Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. ● Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. See the back cover for addresses. Note: Disposal of explanted devices or leads is subject to local, state, and federal regulations. Waste class – B (in accordance with Russian classification) after contact with the body, class A – for waste not on contact with the body.
13. Information about manufacturer of medical product (name, location address, mailing address, phone number, address of place of production)	Organizational and legal form and full name of the legal entity: Medtronic Inc., USA Address (location) of the legal entity: 710 Medtronic Parkway N. E. Minneapolis, MN 55432, USA Phone numbers: Tel.: +1 763 514 4000; Fax: +1 763 514 4879; The e-mail address of the legal entity (if available): http://www.medtronic.com/contact-us/index.htm Places of manufacture of medical device: 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co., Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

14. Classification of medical product: class	Hazard class, depending on the possible risk associated with the use of the medical device - CLASS 3 / III Class of electrical safety, safety when working with lasers and other types of hazardous factors (if any) The software security class: C (according to IEC 62304) Type of contact with the body - permanent contact with soft tissues and interstitial fluids. Type of working part: similar to CF; Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A Class of medical waste: class B - after contact with the body, class A – for waste not on contact with the body.
15. Type of contact with the body	Permanent contact with soft tissues and interstitial fluids .
16. Methods and apparatus for disinfection and sterilization cleaning	Not applicable. Medical device is single use only.
17. Shelf life period	Shelf life -18 months
18. Operating conditions:	Used in medical and preventive treatment facility. Should be used only by or under the supervision of physicians trained. internal body environment 32°C-42°C, humidity 80-100% Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 22°C and 45°C. Sensitivity at nominal conditions as measured at 37°C can vary as much as $\pm 1\%$ per °C from 22°C to 45°C.
19. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. There is no other storage requirements outside of temperature of our IPGs: Store and transport the package between -18 °C and +55 °C.
20. Assembly and installation requirements	Not applicable
21. Warranty	Guaranteed shelf life-18 months. Please see Limited Warranty for information.

22. The list of international regulatory documents/standards, which matches medical product	REFERENCED STANDARDS & REVISIONS	STANDARD TITLE	COMPLIANCE LEVEL
	EN 556-1	Sterilization of medical devices - Requirements for medical devices to be labeled "Sterile"- Part 1: Requirements for terminally sterilized medical devices	Full Compliance
	ISO 15223-1	Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied — Part 1: General requirements Second Edition	Full Compliance
	EN 1041	Information supplied by the manufacturer with medical devices	Full Compliance
	ISO 5841-3	Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
	EN ISO 10993-1	Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing	Full Compliance
	EN ISO 10993-7	Biological evaluation of medical devices - Part 7: Ethylene Oxide sterilization residuals	Full Compliance
	ISO 11135	Sterilization of medical devices – Validation and routine control of Ethylene Oxide sterilization	Full Compliance
	ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2006)	Full Compliance

	ISO 11607-2	Packaging for terminally sterilized medical devices – Part 2: Validation requirements for forming, sealing and assembly processes	Full Compliance
	EN ISO 14971	Medical devices - Application of risk management to medical devices	Full Compliance
	EN 45502-1	Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
	EN 45502-2-1	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia	Full Compliance
	EN 60601-1	General requirements for safety of Medical Electrical devices	Full Compliance (only clause 14 is applicable)
	EN 60601-1-6	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance
	EN 62304	Medical device software. Software life-cycle processes	Full Compliance
	EN 62366	Medical devices - Application of usability engineering to medical devices	Full Compliance
23. Service	Customer Service	Service Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic Inc., USA	

Information for patients:

Toll-free +1 800 551 5544

(7 a.m. – 6 p.m., Monday–Friday, central time)

Fax +1 763 514 1855

Representative company of Medtronic in Moscow

Medtronic LLC,

Russia, 123317 Moscow,

Presnenskaya Naberezhnaya, 10

Tel.: +7 495 580 7377 +7 495 580 7377 FREE

Fax: + 7 495 580 7378

E-mail: info.russia@medtronic.ru

Internet: www.medtronic.eu, www.medtronic.ru

To find the addresses of offices worldwide Medtronic, please go to the following link:

<http://www.medtronic.com>

Appendix #1 to the Technical information

Device parameters

Pacemaker Implantable single chamber Azure XT SR MRI SureScan, model W2SR01

5.1 Emergency settings

Table 11. Emergency VVI settings

Parameter	Selectable values
Pacing Mode	VVI
Lower Rate	70 min ⁻¹
RV Amplitude ^a	6 V
RV Pulse Width ^a	1.5 ms
RV Pace Polarity	Unipolar
V. Blank Post VP	240 ms
Rate Hysteresis	Off
V. Rate Stabilization	Off
MRI SureScan	Off

^a If the programmed RV Amplitude is 8 V, VVI pacing is delivered at 8 V with a pulse width of 1.2 ms.

5.2 Magnet application

When a magnet is placed near the device, the pacing mode changes from the programmed mode to VOO, and the pacing rate changes to 100 min⁻¹ for 5 beats and then changes to 85 min⁻¹ (± 2 min⁻¹) or 65 min⁻¹ (± 2 min⁻¹), as described at the end of this section. Placing a magnet near the device suspends tachyarrhythmia detection. When the magnet is removed, the device returns to its programmed operation. The pacing rate will be 85 min⁻¹ (700 ms) if the device conditions are normal and it will be 65 min⁻¹ (920 ms) if a Recommended Replacement Time (RRT) indicator or a device reset has occurred.

Note: Magnet operation does not occur if telemetry between the device and programmer is established or if the MRI SureScan mode is programmed to On.

5.3 Tachyarrhythmia detection parameters

Table 12. Tachyarrhythmia detection parameters

Parameter	Programmable values	Shipped	Reset
VT Monitor	Monitor  ; Off	Monitor	Off
VT Monitor Interval (Rate) ^a	280; 290 ... 360  ... 500 ms	360 ms	360 ms
RV Sensitivity ^{b,c}	0.45; 0.60 mV ($\pm 50\%$); 0.90; 1.20; 2.00; 2.80; 4.00; 5.60; 8.00; 11.30 mV ($\pm 30\%$) Bipolar: 0.90  mV Unipolar: 2.80  mV	0.90 mV	2.80 mV

^a The measured intervals are truncated to a 10 ms multiple (for example, 457 ms becomes 450 ms). The device uses this truncated interval value when applying the programmed criteria and calculating interval averages.

^b This setting applies to all sensing in this chamber for both tachyarrhythmia detection and bradycardia pacing operations.

^c The device complies with the requirements of ISO 14708-2 when the sensitivity threshold is programmed to 2.0 mV or higher.

5.4 Pacing parameters

Table 13. Modes, rates, and intervals

Parameter	Programmable values	Shipped	Reset
Mode	VVIR [⚙] ; VVI; VOO; OVO	VVI	VVI
Lower Rate ^a	30; 35 ... 60 [⚙] ; 70; 75 ... 150 min ⁻¹ (±2 min ⁻¹)	60 min ⁻¹ (1000 ms)	65 min ⁻¹ (923 ms)

^a The corresponding Lower Rate interval can be calculated as follows: Lower Rate interval (ms) = 60,000/Lower Rate.

Table 14. RV parameters

Parameter	Programmable values	Shipped	Reset
RV Amplitude	0.5; 0.75 ... 1.25 V (+0.125 V / -33%) 1.50 ... 3.5 [⚙] ... 5; 5.5; 6; 8 V (+15% / -33%) ^a	3.5 V	6 V
RV Pulse Width	0.03; 0.06 ms (±10 µs); 0.1; 0.2; 0.3; 0.4 [⚙] ... 1.5 ms (±25 µs)	0.4 ms	1.5 ms
RV Sensitivity ^b	0.45; 0.60 mV (±50%); 0.90 [⚙] ; 1.20; 2.00; 2.80; 4.00; 5.60; 8.00; 11.30 mV (±30%) Unipolar: 2.80 [⚙] mV Bipolar: 0.90 [⚙] mV	0.90 mV	2.80 mV
RV Pace Polarity	Bipolar; Unipolar	Configure ^c	Unipolar
RV Sense Polarity	Bipolar; Unipolar	Configure ^c	Unipolar
RV Lead Monitor	Monitor Only; Adaptive	Monitor Only	Monitor Only
Min Limit	200 [⚙] ; 300; 400; 500 Ω	200 Ω	200 Ω
Max Limit	1000; 1500; 2000; 3000 [⚙] Ω	3000 Ω	3000 Ω

^a When RV Amplitude is 8 V, RV Pulse Width must be less than 1.3 ms.

^b This setting applies to all sensing in this chamber for both tachyarrhythmia detection and bradycardia pacing operations.

^c "Configure" is displayed when the device is automatically configuring the lead polarity at implant. It is not a selectable value.

Table 15. RV Capture Management parameters

Parameter	Programmable values	Shipped	Reset
RV Capture Management	Adaptive [⚙] ; Monitor; Off	Adaptive	Off
RV Amplitude Safety Margin	1.5x; 2.0x [⚙] ; 2.5x; 3.0x	2.0x	2.0x
RV Minimum Adapted Amplitude	1.0; 1.5; 2.0 [⚙] ; 2.5; 3.0; 3.5 V	2 V	2 V
RV Acute Phase Remaining	Off; 30; 60; 90; 120 [⚙] ; 150 days	120 days	120 days

Table 16. Blanking periods

Parameter	Programmable values	Shipped	Reset
-----------	---------------------	---------	-------

V. Blank Post VP	150; 160 ... 200  ... 320 ms (+5; -30 ms)	200 ms	240 ms
V. Blank Post VS	120  ; 130 ... 170; 200; 220; 250; 280; 300; 320 ms (+2; -30 ms)	120 ms	120 ms

Table 17. Rate Response Pacing parameters

Parameter	Programmable values	Shipped	Reset
Upper Sensor Rate	80; 85 ... 130  ... 175 min ⁻¹ (±2 min ⁻¹)	130 min ⁻¹	120 min ⁻¹
ADL Rate	60; 65 ... 95  ... 170 min ⁻¹ (±2 min ⁻¹)	95 min ⁻¹	95 min ⁻¹
Rate Profile Optimization	On  ; Off	On	Off
ADL Response	1; 2; 3  ; 4; 5	3	3
Exertion Response	1; 2; 3  ; 4; 5	3	3
Activity Threshold	Low; Medium Low  ; Medium High; High	Medium Low	Medium Low
Activity Acceleration	15; 30  ; 60 s	30 s	30 s
Activity Deceleration	Exercise  ; 2.5; 5; 10 min	Exercise	5 min
ADL Setpoint	5; 6 ... 40; 42 ... 80	18	18
UR Setpoint	15; 16 ... 40; 42 ... 80; 85 ... 180	40	40

Table 18. Conducted AF Response parameters

Parameter	Programmable values	Shipped	Reset
Conducted AF Response	On; Off 	Off	Off
Response Level	Low; Medium  ; High	Medium	Medium
Maximum Rate	80; 85 ... 110  ... 130 min ⁻¹	110 min ⁻¹	110 min ⁻¹

Table 19. Ventricular Rate Stabilization parameters

Parameter	Programmable values	Shipped	Reset
V. Rate Stabilization	On; Off 	Off	Off
Maximum Rate	80; 85 ... 100  ... 120 min ⁻¹	100 min ⁻¹	120 min ⁻¹
Interval Increment	100; 110 ... 150  ... 400 ms	150 ms	150 ms

Table 20. Sleep parameters

Parameter	Programmable values	Shipped	Reset
Sleep	On; Off 	Off	Off
Sleep Rate	30; 35 ... 50  ; 55; 60; 70; 75 ... 100 min ⁻¹	50 min ⁻¹	50 min ⁻¹
Bed Time	00:00; 00:10 ... 22:00  ... 23:50	22:00	22:00
Wake Time	00:00; 00:10 ... 07:00  ... 23:50	07:00	07:00

Table 21. MRI SureScan parameters

Parameter	Programmable values	Shipped	Reset
MRI SureScan	On; Off	Off	Off

MRI Pacing Mode	VOO; OVO	—	—
MRI Pacing Rate	60; 70; 75; 80 ... 120 min ⁻¹	—	—

Table 22. Additional pacing features

Parameter	Programmable values	Shipped	Reset
Rate Hysteresis	Off [⚡] ; 30; 40 ... 80 min ⁻¹	Off	Off

5.5 Data collection parameters

Table 23. Data collection parameters

Parameter	Programmable values	Shipped	Reset
EGM 1 Source	Can to RVring; RVtip to RVring [⚡] ; RVtip to Can	RVtip to RVring	RVtip to RVring
EGM 1 Range	±1; ±2; ±4; ±8 [⚡] ; ±12; ±16; ±32 mV	±8 mV	±8 mV
EGM 2 Source	Can to RVring; RVtip to RVring; RVtip to Can [⚡]	RVtip to Can	RVtip to Can
EGM 2 Range	±1; ±2; ±4; ±8 [⚡] ; ±12; ±16; ±32 mV	±8 mV	±8 mV
EGM 3 Source	RVtip to RVring; Can to RVring [⚡] ; RVtip to Can	Can to RVring	Can to RVring
EGM 3 Range	±1; ±2; ±4; ±8 [⚡] ; ±12; ±16; ±32 mV	±8 mV	±8 mV
Monitored	EGM1 and EGM2 [⚡] ; EGM1 and EGM3; EGM2 and EGM3	EGM1 and EGM2	EGM1 and EGM2
Pre-arrhythmia EGM	Off [⚡] ; On - 1 month; On - 3 months; On Continuous	Off	Off
OptiVol Threshold ^{a,b}	30; 40; 50; 60 [⚡] ... 160; 170; 180	60	60
Device Date/Time ^b	(select Time Zone)	—	—
Holter Telemetry	Off [⚡] ; 0.5; 1; 2; 4; 8; 16; 24; 36; 46 hr	Off	Off
Wireless Telemetry with Monitor	On [⚡] ; Off	On	On ^d

^a The OptiVol Threshold settings are available in the Data Collection Setup screen only when the Wireless Telemetry with Monitor parameter is programmed to Off.

^b Decreasing the OptiVol Threshold will make the device more sensitive to changes in the patient's thoracic fluid status. Increasing the OptiVol Threshold could delay or prevent device observation of significant changes in the patient's thoracic fluid status.

^c The times and dates stored in episode records and other data are determined by the Device Date/Time clock.

^d The reset value may be set to Off if there is an issue with wireless communication that requires it to be disabled.

5.6 Medtronic CareAlert parameters

Table 24. Clinical Management Alerts

Parameter	Programmable values	Shipped	Reset
Monitored VT Episode Detected	Off [⚡] ; On	Off	Off
OptiVol 2.0 Fluid Settings...			
OptiVol Alert Enable	Off [⚡] ; On	Off	Off

OptiVol Threshold ^a	30; 40; 50; 60  ... 100; 120 ... 180	60	60
--------------------------------	---	----	----

^a Decreasing the OptiVol Threshold makes the device more sensitive to changes in the patient's thoracic fluid status. Increasing the OptiVol Threshold could delay or prevent device observation of significant changes in the patient's thoracic fluid status.

Table 25. Lead/Device Integrity Alerts

Parameter	Programmable values	Shipped	Reset
Low Battery Voltage RRT	On  ; Off	On	Off
Lead Impedance Out of Range...			
Lead Impedance			
RV Pacing Enable	On  ; Off	On	Off
RV Pacing Less than	200  ; 300; 400; 500 Ω	200 Ω	200 Ω
RV Pacing Greater than	1000; 1500; 2000; 3000  Ω	3000 Ω	3000 Ω
Capture Management High Threshold...			
High Threshold			
RV Capture Enable ^a	Off  ; On	Off	Off

^a If programmed to On, alert notification is sent if RV capture management has measured high thresholds for 3 consecutive days.

5.7 System test parameters

Table 26. System test parameters

Parameter	Selectable values
Pacing Threshold Test parameters	
Test Type	Amplitude; Pulse Width
Chamber	RV
Decrement after	2; 3 ... 15 pulses
RV Pace Polarity	Unipolar; Bipolar
Mode ^a	VVI; VOO
Lower Rate	30; 35 ... 60; 70; 75 ... 150 min ⁻¹
RV Amplitude	0.25; 0.5 ... 5; 5.5; 6; 8 V
RV Pulse Width	0.03; 0.06; 0.1; 0.2 ... 1.5 ms
V. Pace Blanking	150; 160 ... 320 ms
Sensing Test parameters	
Mode ^a	VVI; OVO
Lower Rate	30; 35 ... 60; 70; 75 ... 120 min ⁻¹

^a The selectable values for this parameter depend on the programmed pacing mode.

5.8 EP study parameters

Table 27. Fixed Burst induction parameters

Parameter	Selectable values
Interval	100; 110 ... 600  ms

Amplitude	1; 2; 3; 4 ; 5; 6; 8 V
Pulse Width	0.10; 0.20 ... 0.50 ... 1.50 ms

Table 28. PES induction parameters

Parameter	Selectable values
#S1	1; 2 ... 8 ... 15
S1S1	100; 110 ... 600 ... 2000 ms
S1S2	Off; 100; 110 ... 400 ... 600 ms
S2S3	Off ; 100; 110 ... 600 ms
S3S4	Off ; 100; 110 ... 600 ms
Amplitude	1; 2; 3; 4 ; 5; 6; 8 V
Pulse Width	0.10; 0.20 ... 0.50 ... 1.50 ms

Table 29. Shared manual ATP therapy parameters

Parameter	Selectable values
Minimum Interval (ventricular ATP)	150; 160 ... 200 ... 400 ms
Amplitude	1; 2 ... 6 ; 8 V
Pulse Width	0.10; 0.20 ... 1.50 ms

Table 30. Manual Ramp therapy parameters

Parameter	Selectable values
Chamber	RV
RV Ramp therapy parameters	
# Pulses	1; 2 ... 6 ... 15
%RR Interval	50; 53; 56; 59; 63; 66 ... 84; 88; 91; 94; 97 %
Dec/Pulse	0; 10 ; 20; 30; 40 ms

Table 31. Manual Burst therapy parameters

Parameter	Selectable values
# Pulses	1; 2 ... 8 ... 15
%RR Interval	50; 53; 56; 59; 63; 66 ... 84; 88 ; 91; 94; 97 %

Table 32. Manual Ramp+ therapy parameters

Parameter	Selectable values
# Pulses	1; 2; 3 ... 15
R-S1 (%RR)	50; 53; 56; 59; 63; 66 ... 75 ... 84; 88; 91; 94; 97%
S1-S2 (%RR)	50; 53; 56; 59; 63; 66; 69 ... 84; 88; 91; 94; 97%
S2-SN (%RR)	50; 53; 56; 59; 63; 66 ... 84; 88; 91; 94; 97%

5.9 Nonprogrammable parameters

Table 33. Nonprogrammable parameters

Parameter	Value
Premature event threshold for counting PVCs and Runs of PVCs	95%
Hardware parameters	
Pacing rate limit (protective feature)	200 min ⁻¹
Input impedance	150 kΩ minimum
Effective pacing capacitance	4 μF (±10%)
Recommended Replacement Time (RRT)	
Battery Voltage Threshold	≤ 2.63 V

Table 34. Nonprogrammable parameters for the MRI SureScan mode

Parameter	Value
Pacing amplitude	Programmed pacing amplitude value when >5 V; 5 V when programmed pacing amplitude value is ≤5 V
Pulse width	Programmed pulse width value when >1 ms; 1 ms when programmed pulse width value is ≤1 ms
Sensitivity	Programmed value
Input impedance	150 kΩ
Pacing rate limit	200 min ⁻¹
Effective pacing capacitance	4 μF (±10%)
Blanking period	
OVO mode	Programmed blanking period value
VOO mode	—

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
Addendum to Device Manual

1. Name of medical product	Pacemaker implantable single-chamber Azure S SR MRI SureScan, model W3SR01
2.1. Prescription of a medical product	Medical device is a multiprogrammable cardiac device, which monitors and regulates the patient's heart rate, by performing a one-chamber rate-responsive adaptive stimulation with bradycardia.
2.2 Description of the principles underlying the work of the medical products and their features	Feature summary Pacing features Auto-adjusting sensitivity – This feature automatically adjusts the sensitivity thresholds after specific paced events and sensed events occur. Automatic polarity configuration – This device uses lead impedance measurements to automatically configure pacing and sensing polarities during Implant Detection. MRI SureScan – This feature allows patients with an implanted MRI SureScan system, including the device and lead, to have a safe MRI procedure if the requirements provided in the MRI technical manual are followed. Rate Profile Optimization – The goal of Rate Profile Optimization is to ensure that the rate response remains appropriate for the full range of patient activities. This feature monitors the patient's daily and monthly sensor rate profiles and adjusts the rate response curves over time to achieve a prescribed target rate profile.

Cindi A. Bigelbach
12/10/2019

Rate-responsive pacing – This feature varies the pacing rate in response to the patient’s physical motion as detected by the activity sensor of the device.

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Inductive telemetry communication system – The Medtronic inductive telemetry communication system is used with the clinician programmer to interrogate and program the device. It can also be used to interrogate the device for

remote monitoring, if available. This system uses short-range communication that protects patient information and device data.

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Position:	I	II	III	IV	V
Category:	Chamber(s) Paced	Chamber(s) Sensed	Response to Sensing	Rate Modula- tion	Multisite Pac- ing ^a
	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None T = Triggered I = Inhibited D = Dual (T + I)	O = None R = Rate mod- ulation	O = None A = Atrium V = Ventricle D = Dual (A + V)
Manufacturers' designation only:	S = Single ^b (A or V)	S = Single ^b (A or V)			

^a Medtronic devices do not use the Multisite Pacing code.

^b The programmer displays A or V (not S) for chambers paced and sensed.

3. Conditions for use

In medical and preventive treatment facility.

The pacemaker should be used only by or under the supervision of physicians trained.

4. Composition and list of parts and accessories*





Implantable pacemaker single-chamber Azure S SR MRI SureScan, model W3SR01 consisting of:

- implantable pacemaker single-chamber Azure S SR MRI SureScan, model W3SR01;

- torque wrench

-Product documentation includes - Device manual; Limited Warranty; Medical procedures and EMI Warnings and Precautions manual; MRI Technical manual ; Explanation of the symbols IPG and CRT-P devices; Stickers; Registration card; Medtronic Web Card; CD-ROM.


	* Appendix to the Device Manual isn't included in package content, provided to customer separately on the territory of Russian Federation, SW030 Programmer software guide and Reference Manual isn't included in package content, visit http://manuals.medtronic.com/						
5. Technical details of medical product	Pacemaker implantable single-chamber Azure S SR MRI SureScan, model W3SR01:						
	Volume ^a	12.25 cm ³					
	Mass	22.5 g					
	H x W x D ^b	42.6 mm x 50.8 mm x 7.4 mm					
	Radiopaque ID ^c	RNA					
	Medtronic identifier						
	Surface area of titanium device can	33.48 cm ²					
	Materials in contact with human tissue ^d	Titanium, polyurethane, silicone rubber					
	Battery	Lithium-hybrid CFx silver vanadium oxide					
	^a Volume with connector holes unplugged.						
^b Grommets may protrude slightly beyond the can surface.							
^c The radiopaque ID, which includes a Medtronic-identifier symbol, can be viewed in a fluoroscopic image of the device.							
^d These materials have been successfully tested for the ability to avoid biological incompatibility. The device does not produce an injurious temperature in the surrounding tissue during normal operation.							
		Height (mm)	Width (mm)	Case thickness (mm)	Connector Thickness (mm)	Mass (g)	Volume (cm ³)
Single Chamber Mode		42.6	50.8	7.4	9.9	22.5	12.25
Tolerance		+/-10%	+/-10%	+/-10%	+/-10%	+/-10%	+/-10%
Radiopaq ID drawing and location:							

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
	The Astra and Azure MRI Systems may be comprised of the following components: • CareLink Programmer Model 2090 • CareLink Encore Programmer Model 29901 • Conexus Activator Model 27901 • MyCareLink Monitor Model 24952 • MyCarelink Smart, Model 25000 (Astra only) • 2020B CareLink Express Monitor • MyCareLink Heart Model 31302 and Model 31301 (Azure only) • Commercially available MRI SureScan Right Atrial Pacing Leads¹ • Commercially available MRI SureScan Left Ventricular Pacing Leads² • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected) • The same commercially available implant support instruments and accessories used with the predecessor Advisa DR device systems. <i>Commercially Available Accessories</i> The Azure/Astra system will be compatible with the following commercially available system components: • Model 2290 pacing system analyzer • Model 9466 Patient Magnet • Model DR220 Holter Monitor • Model PA96000 Patient Assistant
--	--

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


		accordance with IEC 61000-4-6	LPM of HF devices		
		The emitted radio-frequency electromagnetic field according to IEC 61000-4-3	10 V (RMS) in the band from 150 kHz to 80 MHz at the frequencies allocated for the LPM of high-frequency devices	2.3 V (rms value)	1 V (rms value)
				13 V/m	10 V/m
		Noise immunity			
		The test for Noise immunity	Test Level according to IEC 60601	Level Conformity	Electromagnetic environment - instructions
		Electrostatic Discharge (ESR) according to IEC 61000-4-2	± 6 kV – contact discharge ± 8 kV – air discharge	± 6 kV – contact discharge ± 8 kV – air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%.
		A magnetic field Industrial Frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of the industrial frequency must correspond to the typical conditions of a commercial or hospital environment.
		Radio frequency Electromagnetic field IEC 61000-4-3	3 V / m In the band from 80 MHz to 2.5 GHz	3 V/m	However, the readings of the weights built into the bed do not have a critical value for the operation of a medical device: D = 1,2√P in the band from 26 to 800 MHz D = 2.3√P in the band from 800 MHz to 2.5 GHz,

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
	NOTE 1 Higher frequency bands are applied at 80 MHz and 800 MHz. NOTE 2: This manual is not applicable in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.
8. Materials used for production of medical device	Implantable pacemaker Titanium grade 1, Silicone rubber, Polyurethane Torque wrench Polyetherimide, Corrosion resistant steel
9. Design and appearance, photo or image	Pacemaker implantable single-chamber Azure S SR MRI SureScan, model W3SR01 consisting of: Pacemaker implantable single-chamber Azure S SR MRI SureScan, model W3SR01  Torque wrench

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



Pacemaker implantable single-chamber Azure S SR MRI SureScan, model W3SR01 (shelf box and single package labels)

Azure™ S SR MRI
SureScan™

Model W3SR01
SN WWW000000W
2000-00-00
YYYY-MM-DD
2000-00-00
YYYY-MM-DD
REF W3SR01


VVI

60 min⁻¹

0.9 mV

1 3.5 V ♂
-0.4 ms ♂

AUTO

AUTO

85 min⁻¹

STERILE EO

-18 °C
0 °F

55 °C
131 °F

www.medtronic.com/manuals

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands

Manufactured in: Tolochenaz, Switzerland

Product of Switzerland
www.medtronic.com

CE0123/1856

8970311A001 B - M967976A001 A


Azure™ S SR MRI
SureScan™


Model W3SR01
SN WWW000000W
2000-00-00
YYYY-MM-DD
2000-00-00
YYYY-MM-DD
REF W3SR01


VVI

60 min⁻¹

0.9 mV


1 3.5 V ♂
-0.4 ms ♂

AUTO

AUTO

85 min⁻¹

FCC ID: LFSBLEIMPLANT

IC: 3408D-BLEIMPLANT

STERILE EO

-18 °C
0 °F

55 °C
131 °F

www.medtronic.com/manuals

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA

Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands

Manufactured in: Tolochenaz, Switzerland

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
12. Requirements for the protection of the environment in the application of medical products

Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See <http://recycling.medtronic.com> for instructions on proper disposal of this product.
Consider the following information related to device explant and disposal:

	● To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. ● Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. ● Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. ● Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. See the back cover for addresses. Note: Disposal of explanted devices or leads is subject to local, state, and federal regulations. Waste class – B (in accordance with Russian classification) after contact with the body, class A – for waste not on contact with the body.
13. Information about manufacturer of medical product (name, location address, mailing address, phone number, address of place of production)	Organizational and legal form and full name of the legal entity: Medtronic Inc., USA Address (location) of the legal entity: 710 Medtronic Parkway N. E. Minneapolis, MN 55432, USA Phone numbers: Tel.: +1 763 514 4000; Fax: +1 763 514 4879; The e-mail address of the legal entity (if available): http://www.medtronic.com/contact-us/index.htm Places of manufacture of medical device: 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co., Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA
14. Classification of medical product: class	Hazard class, depending on the possible risk associated with the use of the medical device - CLASS 3 / III Class of electrical safety, safety when working with lasers and other types of hazardous factors (if any) The software security class: C (according to IEC 62304) Type of contact with the body - permanent contact with soft tissues and interstitial fluids. Type of working part: similar to CF; Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A

	Class of medical waste: class B - after contact with the body, class A – for waste not on contact with the body.		
15. Type of contact with the body	Permanent contact with soft tissues and interstitial fluids.		
16. Methods and apparatus for disinfection and sterilization cleaning	Not applicable. Medical device is single use only.		
17. Shelf life period	Shelf life -18 months		
18. Operating conditions:	Used in medical and preventive treatment facility. Should be used only by or under the supervision of physicians trained. internal body environment 32°C-42°C, humidity 80-100% Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 22°C and 45°C. Sensitivity at nominal conditions as measured at 37°C can vary as much as $\pm 1\%$ per °C from 22°C to 45°C.		
19. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. There is no other storage requirements outside of temperature of our IPGs: Store and transport the package between -18 °C and +55 °C.		
20. Assembly and installation requirements	Not applicable		
21. Warranty	Guaranteed shelf life-18 months. Please see Limited Warranty for information.		
22. The list of international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS	STANDARD TITLE	COMPLIANCE LEVEL

, which matches medical product	EN 556-1	Sterilization of medical devices - Requirements for medical devices to be labeled "Sterile"- Part 1: Requirements for terminally sterilized medical devices	Full Compliance
	ISO 15223-1	Medical devices — Symbols to be used with medical device labels, labeling, and information to be supplied — Part 1: General requirements Second Edition	Full Compliance
	EN 1041	Information supplied by the manufacturer with medical devices	Full Compliance
	ISO 5841-3	Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
	EN ISO 10993-1	Biological Evaluation of Medical Devices - Part 1: Evaluation and Testing	Full Compliance
	EN ISO 10993-7	Biological evaluation of medical devices - Part 7: Ethylene Oxide sterilization residuals	Full Compliance
	ISO 11135	Sterilization of medical devices – Validation and routine control of Ethylene Oxide sterilization	Full Compliance
	ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2006)	Full Compliance
	ISO 11607-2	Packaging for terminally sterilized medical devices – Part 2: Validation	Full Compliance

		requirements for forming, sealing and assembly processes	
	EN ISO 14971	Medical devices - Application of risk management to medical devices	Full Compliance
	EN 45502-1	Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
	EN 45502-2-1	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia	Full Compliance
	EN 60601-1	General requirements for safety of Medical Electrical devices	Full Compliance (only clause 14 is applicable)
	EN 60601-1-6	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance
	EN 62304	Medical device software. Software life-cycle processes	Full Compliance
	EN 62366	Medical devices - Application of usability engineering to medical devices	Full Compliance
23. Service	Customer	Service Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544	

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EGM 1 Range	± 1 ; ± 2 ; ± 4 ; ± 8  ; ± 12 ; ± 16 ; ± 32 mV	± 8 mV	± 8 mV
EGM 2 Source	Can to RVring; RVtip to RVring; RVtip to Can 	RVtip to Can	RVtip to Can
EGM 2 Range	± 1 ; ± 2 ; ± 4 ; ± 8  ; ± 12 ; ± 16 ; ± 32 mV	± 8 mV	± 8 mV
EGM 3 Source	RVtip to RVring; Can to RVring  ; RVTip to Can	Can to RVring	Can to RVring
EGM 3 Range	± 1 ; ± 2 ; ± 4 ; ± 8  ; ± 12 ; ± 16 ; ± 32 mV	± 8 mV	± 8 mV
Monitored	EGM1 and EGM2  ; EGM1 and EGM3; EGM2 and EGM3	EGM1 and EGM2	EGM1 and EGM2
Pre-arrhythmia EGM	Off  ; On - 1 month; On - 3 months; On Continuous	Off	Off
Device Date/Time ^a	(select Time Zone)	—	—
Holter Telemetry	Off  ; 0.5; 1; 2; 4; 8; 16; 24; 36; 46 hr	Off	Off
Wireless Telemetry with Monitor	On  ; Off	On	On ^b

^a The times and dates stored in episode records and other data are determined by the Device Date/Time clock.

^b The reset value may be set to Off if there is an issue with wireless communication that requires it to be disabled.

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
Addendum to Device Manual

1. Name of medical product	Pacemaker implantable single-chamber Astra XT SR MRI SureScan, model X2SR01
2.1. Prescription of a medical product	Medical device is a multiprogrammable cardiac device, which monitors and regulates the patient's heart rate, by performing a one-chamber rate-responsive adaptive stimulation with bradycardia.
2.2 Description of the principles underlying the work of the medical products and their features	Feature summary Pacing features Auto-adjusting sensitivity – This feature automatically adjusts the sensitivity thresholds after specific paced events and sensed events occur. Automatic polarity configuration – This device uses lead impedance measurements to automatically configure pacing and sensing polarities during Implant Detection. Conducted AF Response – This feature regularizes the ventricular rhythm during conducted AT/AF by modifying the pacing rate on a beat-by-beat basis to closely match the patient's average ventricular response. MRI SureScan – This feature allows patients with an implanted MRI SureScan system, including the device and lead, to have a safe MRI procedure if the requirements provided in the MRI technical manual are followed.

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	Position:	I	II	III	IV	V
	Category:	Chamber(s) Paced	Chamber(s) Sensed	Response to Sensing	Rate Modula- tion	Multisite Pac- ing ^a
		O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None T = Triggered I = Inhibited D = Dual (T + I)	O = None R = Rate mod- ulation	O = None A = Atrium V = Ventricle D = Dual (A + V)
	Manufacturers' designation only:	S = Single ^b (A or V)	S = Single ^b (A or V)			
^a Medtronic devices do not use the Multisite Pacing code.						
^b The programmer displays A or V (not S) for chambers paced and sensed.						
3. Conditions for use	In medical and preventive treatment facility. The pacemaker should be used only by or under the supervision of physicians trained.					
4. Composition and list of parts and accessories*	Implantable pacemaker single-chamber Astra XT SR MRI SureScan, model X2SR01 consisting of: - implantable pacemaker single-chamber Astra XT SR MRI SureScan, model X2SR01; - torque wrench -Product documentation includes - Device manual; Limited Warranty; Medical procedures and EMI Warnings and Precautions manual; MRI Technical manual ; Explanation of the symbols IPG and CRT-P devices; Stickers; Registration card; Medtronic Web Card; CD-ROM. * Appendix to the Device Manual isn't included in package content, provided to customer separately on the territory of Russian Federation, SW030 Programmer software guide and Reference Manual isn't included in package content, visit http://manuals.medtronic.com/					
5. Technical details of medical product	Pacemaker implantable single-chamber Astra XT SR MRI SureScan, model X2SR01:					
	Volume ^a		12.25 cm ³			
	Mass		22.5 g			
	H x W x D ^b		42.6 mm x 50.8 mm x 7.4 mm			
	Radiopaque ID ^c		RNA			

Medtronic identifier

Surface area of titanium device can

33.48 cm²

Materials in contact with human tissue^d

Titanium, polyurethane, silicone rubber

Battery

Lithium-hybrid CFx silver vanadium oxide

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




Radiopaq ID drawing and location:

SECTION E-E

All measurements and tolerances for radiopaque drawing is in inches.
Unless otherwise specified, default tolerance is +/- 0.006 inches

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Only the physician programmer, using inductive telemetry communication system is capable of programming the Astra and Azure IPGs. Patients may use the optional Medtronic patient monitor to gather information from their implanted devices and to send information to their physicians through the Medtronic CareLink Network (patient-initiated transmissions). The alternative, optional, wireless Bluetooth Low Energy (BTLE) mobile application (not subject to medical device regulation) and BTLE enabled Medtronic patient monitor (available for Azure models only) allow device initiated (CareAlerts) and remote scheduled transmissions without patient initiation. The IPG will only respond to a restricted set of commands, such as interrogation commands. These commands do not allow therapy parameters to be modified. Data sent from the implanted device through the wireless communication system will be encrypted in the IPG and will pass through the mobile app to the CareLink server where it is decrypted. The transmissions between the IPG and CareLink server are encrypted end-to-end.

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
	Lead Model	Brand Name	
	5086MRI	CapSureFix MRI SureScan Pace/Sense Lead	
	5076	CapSureFix Novus MRI SureScan Pace/Sense Lead	
	4076	CapSureFix Novus MRI SureScan Lead	
	4074/4574	CapSure Sense MRI SureScan Lead	
	3830	SelectSecure Model Pacing Lead	
	5554	CapSure Z Novus Model 5554 pacing lead	
	5054	CapSure Z Novus Model 5054 pacing lead	
7. EMC Information	Electromagnetic emission		
	Medical device are intended for use in the electromagnetic environment defined below. Buyer or user medical device should ensure its use in the specified environment.		
	Electromagnetic emission test	Conformity	Electromagnetic environment - instructions
	Radio interference according to CISPR 11	Group 1	Medical device should emit electromagnetic energy to perform the main function.
	Radio interference to CISPR 11	Class B	Medical device are suitable for use in all locations.
	Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use	

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

		accordance with IEC 61000-4-6	LPM of HF devices		
		The emitted radio-frequency electromagnetic field according to IEC 61000-4-3	10 V (RMS) in the band from 150 kHz to 80 MHz at the frequencies allocated for the LPM of high-frequency devices	2.3 V (rms value)	1 V (rms value)
				13 V/m	10 V/m
	Noise immunity				
	The test for Noise immunity	Test Level according to IEC 60601	Level Conformity	Electromagnetic environment - instructions	
	Electrostatic Discharge (ESR) according to IEC 61000-4-2	± 6 kV – contact discharge ± 8 kV – air discharge	± 6 kV – contact discharge ± 8 kV – air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%.	
	A magnetic field Industrial Frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of the industrial frequency must correspond to the typical conditions of a commercial or hospital environment.	
	Radio frequency Electromagnetic field IEC 61000-4-3	3 V / m In the band from 80 MHz to 2.5 GHz	3 V/m	However, the readings of the weights built into the bed do not have a critical value for the operation of a medical device: D = 1,2√P in the band from 26 to 800 MHz D = 2.3√P in the band from 800 MHz to 2.5 GHz,	

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
	NOTE 1 Higher frequency bands are applied at 80 MHz and 800 MHz. NOTE 2: This manual is not applicable in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.
8. Materials used for production of medical device	Implantable pacemaker Titanium grade 1, Silicone rubber, Polyurethane Torque wrench Polyetherimide, Corrosion resistant steel
9. Design and appearance, photo or image	Pacemaker implantable single-chamber Astra XT SR MRI SureScan, model X2SR01 consisting of: Pacemaker implantable single-chamber Astra XT SR MRI SureScan, model X2SR01  Torque wrench


	
10. Information on labelling, packaging	Sterile barrier systems constructed of uncoated PETG (Glidex) trays, sealed to Oliver 18B coated Tyvek® 1073B porous lid stock and sterilized with ethylene oxide gas. This sterile barrier material configuration is common to many CRDM device models across different product families. Much of the testing that supports the use of packaging materials and sterile package assemblies is specific to the packaging materials and assemblies, not specific products. Product specific evidence for sterilization qualification, package performance qualification, labeling, and product specific requirements can be found in the DHF for each individual product. CRDM Sterile Barrier systems constructed of uncoated PETG trays, sealed to Oliver 18B coated Tyvek® 1073B porous lid stock and sterilized with ethylene oxide gas include the IPG/ICD/CRT Device Package and the Dual Entry Leads Package. Images of these packages follow:



Package dimensions:
1. Sterile package lid (in inches)
Inner package lid:


L

Outer package lid:

Sterile Package dimensions:**Outer package Drawing with Dimensions:**

2

	*Tolerances of dimensions for each package are $\pm 5\%$, unless otherwise specified Labeling of medical devices packages: Pacemaker implantable single-chamber Astra XT SR MRI SureScan, model X2SR01 (shelf box and single package labels)
--	---

Medtronic

Astra™ XT SR MRI
SureScan™

VVIR IS-1

MR

1x

Adaptive

VVI 60 min⁻¹

0.9 mV

1 3.5 V 1/2
-- 0.4 ms 1/2

AUTO

85 min⁻¹

STERILE EO

56 °C 131 °F

-18 °C 0 °F

www.medtronic.com/manuals

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands
Manufactured in: Tolochenaz, Switzerland
Product of Switzerland
www.medtronic.com

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

11. Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)

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12. Requirements for the protection of the environment in the application of medical products	Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See http://recycling.medtronic.com for instructions on proper disposal of this product. Consider the following information related to device explant and disposal: ● To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. ● Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. ● Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. ● Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. See the back cover for addresses. Note: Disposal of explanted devices or leads is subject to local, state, and federal regulations. Waste class – B (in accordance with Russian classification) after contact with the body, class A – for waste not on contact with the body.
13. Information about manufacturer of medical product (name, location address, mailing address, phone number, address of place of production)	Organizational and legal form and full name of the legal entity: Medtronic Inc., USA Address (location) of the legal entity: 710 Medtronic Parkway N. E. Minneapolis, MN 55432, USA Phone numbers: Tel.: +1 763 514 4000; Fax: +1 763 514 4879; The e-mail address of the legal entity (if available): http://www.medtronic.com/contact-us/index.htm Places of manufacture of medical device: 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co., Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

14. Classification of medical product: class	Hazard class, depending on the possible risk associated with the use of the medical device - CLASS 3 / III Class of electrical safety, safety when working with lasers and other types of hazardous factors (if any) The software security class: C (according to IEC 62304) Type of contact with the body - permanent contact with soft tissues and interstitial fluids. Type of working part: similar to CF; Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A Class of medical waste: class B - after contact with the body, class A – for waste not on contact with the body.
15. Type of contact with the body	Permanent contact with soft tissues and interstitial fluids
16. Methods and apparatus for disinfection and sterilization cleaning	Not applicable. Medical device is single use only.
17. Shelf life period	Shelf life -18 months
18. Operating conditions:	Used in medical and preventive treatment facility. Should be used only by or under the supervision of physicians trained. internal body environment 32°C-42°C, humidity 80-100% Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 22°C and 45°C. Sensitivity at nominal conditions as measured at 37°C can vary as much as $\pm 1\%$ per °C from 22°C to 45°C.
19. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. There is no other storage requirements outside of temperature of our IPGs: Store and transport the package between -18 °C and +55 °C.
20. Assembly and installation requirements	Not applicable
21. Warranty	Guaranteed shelf life-18 months. Please see Limited Warranty for information.

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Conducted AF Response	On; Off	Off	Off
Response Level	Low; Medium ; High	Medium	Medium
Maximum Rate	80; 85 ... 110 ... 130 min ⁻¹	110 min ⁻¹	110 min ⁻¹

Table 19. Ventricular Rate Stabilization parameters

Parameter	Programmable values	Shipped	Reset
V. Rate Stabilization	On; Off	Off	Off
Maximum Rate	80; 85 ... 100 ... 120 min ⁻¹	100 min ⁻¹	120 min ⁻¹
Interval Increment	100; 110 ... 150 ... 400 ms	150 ms	150 ms

Table 20. Sleep parameters

Parameter	Programmable values	Shipped	Reset
Sleep	On; Off	Off	Off
Sleep Rate	30; 35 ... 50 ; 55; 60; 70; 75 ... 100 min ⁻¹	50 min ⁻¹	50 min ⁻¹
Bed Time	00:00; 00:10 ... 22:00 ... 23:50	22:00	22:00
Wake Time	00:00; 00:10 ... 07:00 ... 23:50	07:00	07:00

Table 21. MRI SureScan parameters

Parameter	Programmable values	Shipped	Reset
MRI SureScan	On; Off	Off	Off
MRI Pacing Mode	VOO; OVO	—	—
MRI Pacing Rate	60; 70; 75; 80 ... 120 min ⁻¹	—	—

Table 22. Additional pacing features

Parameter	Programmable values	Shipped	Reset
Rate Hysteresis	Off ; 30; 40 ... 80 min ⁻¹	Off	Off

5.5 Data collection parameters

Table 23. Data collection parameters

Parameter	Programmable values	Shipped	Reset
EGM 1 Source	Can to RVring; RVtip to RVring ; RVtip to Can	RVtip to RVring	RVtip to RVring
EGM 1 Range	±1; ±2; ±4; ±8 ; ±12; ±16; ±32 mV	±8 mV	±8 mV
EGM 2 Source	Can to RVring; RVtip to RVring; RVtip to Can	RVtip to Can	RVtip to Can

EGM 2 Range	± 1 ; ± 2 ; ± 4 ; ± 8 ; ± 12 ; ± 16 ; ± 32 mV	± 8 mV	± 8 mV
EGM 3 Source	RVtip to RVring; Can to RVring ; RVTip to Can	Can to RVring	Can to RVring
EGM 3 Range	± 1 ; ± 2 ; ± 4 ; ± 8 ; ± 12 ; ± 16 ; ± 32 mV	± 8 mV	± 8 mV
Monitored	EGM1 and EGM2 ; EGM1 and EGM3; EGM2 and EGM3	EGM1 and EGM2	EGM1 and EGM2
Pre-arrhythmia EGM	Off ; On - 1 month; On - 3 months; On Continuous	Off	Off
OptiVol Threshold ^a	30; 40; 50; 60 ...160; 170; 180 Ω -days	60 Ω -days	60 Ω -days
Device Date/Time ^b	(select Time Zone)	—	—
Holter Telemetry	Off ; 0.5; 1; 2; 4; 8; 16; 24; 36; 46 hr	Off	Off

^a Decreasing the OptiVol Threshold will make the device more sensitive to changes in the patient's thoracic fluid status. Increasing the OptiVol Threshold could delay or prevent device observation of significant changes in the patient's thoracic fluid status.

^b The times and dates stored in episode records and other data are determined by the Device Date/Time clock.

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

CE0123

M970311A001 B - M967982A001 A

Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands
Manufactured in Tolochenaz, Switzerland
Product of Switzerland
www.medtronic.com

Astra™ XT SR MRI
SureScan™

X2SR01

WWW000000W


2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

X2SR01


VVI


60 min⁻¹


AUTO


0.9 mV

85 min⁻¹

3.5 V
0.4 ms

AUTO

FCC ID: LF58BLEIMPLANT
IC: 3408D-BLEIMPLANT

STERILE EO

-18 °C
0 °F

55 °C
131 °F

www.medtronic.com/manuals


Medtronic, Inc., 710 Medtronic Parkway, Minneapolis, MN 55432 USA
Medtronic B.V., Earl Bakkenstraat 10, 6422 PJ Heerlen, The Netherlands
Manufactured in Tolochenaz, Switzerland
www.medtronic.com

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