Yuri K. Carlson

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI Quad CRT-D SureScan, model DTMA2QQ consisting of: – implantable cardioverter defibrillator CLARIA MRI Quad CRT-D SureScan, model DTMA2QQ; – torque wrench; – Device Manual CLARIA MRI Quad CRT-D SureScan DTMA2QQ; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliat 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA

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	Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Claria MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device) ● air embolism⁴ ● bleeding

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10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI Quad CRT-D SureScan, model DTMA2QQ	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper

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The mass of the package with the device, g: 121 $\pm$ 5%;

Secondary Package:


Dimensions, L x W x H: 17.54 cm $\pm$ 0.8 cm x 14.31 cm $\pm$ 0.7 cm x 4.51 cm $\pm$ 0.18 cm;





The mass of the package with the device, g: 513 $\pm$ 5%;

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
Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
	Reorder number
	Serial number
	Use by
	Consult instructions for use
	Date of manufacture
	Package contents

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

and validating software used in a finished medical device	
13. Sterilization method, Sterilization Validation	
14. Procedure and conditions for the	

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	EN 45502-2-2 Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance		
	EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14		
	EN 60601-1-6 Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance		
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance		
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance		
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. <table><tr><td>Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855</td><td>Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com</td></tr></table>		Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855	Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu , www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com
Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855	Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu , www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com			

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

Minori Carlsson
 Minori Carlsson
 Senior Regulatory Affairs Director
 Medtronic, Inc.

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI CRT-D SureScan, model DTMA2D4 consisting of: – implantable cardioverter defibrillator CLARIA MRI CRT-D SureScan, model DTMA2D4; – torque wrench; – Device Manual CLARIA MRI CRT-D SureScan DTMA2D4; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA

	Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Claria MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device) ● air embolism ⁴ ● bleeding

- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)
- tissue necrosis
- valve damage (particularly in fragile hearts)

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
- the Claria MRI device
- a SureScan right atrial pacing lead or a Model 6725 pin plug for the right atrial port
- a SureScan left ventricular pacing lead
- a SureScan defibrillation lead.

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	• For patients whose device will be programmed to an asynchronous pacing mode when MRI SureScan is programmed to On, no diaphragmatic stimulation is present when the paced leads have a pacing output of 5.0 V and a pulse width of 1.0 ms. Note: The LV lead is not paced during SureScan operation so the presence of diaphragmatic stimulation on the LV lead at a pacing output of 5.0 V and a pulse width of 1.0 ms does not need to be considered. Caution: For pacemaker-dependent patients, it is not recommended to perform an MRI scan if the right ventricular (RV) lead pacing capture threshold is greater than 2.0 V at 0.4 ms. A higher pacing capture threshold may indicate an issue with the implanted lead. Patient monitoring and rescue requirements • Continuous patient monitoring is required while MRI SureScan mode is programmed to On. • In the event that patient rescue is required, an external defibrillator must be immediately available. Note: For radiology requirements for an MRI scan, refer to the MRI technical manual. Pre-implant consideration Patient evaluation for the implant of Claria MRI different models should include the following consideration about a concomitant implant with a neurostimulator: Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads; • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected); • The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI CRT-D SureScan, model DTMA2D4	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper

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The mass of the package with the device, g: 120 $\pm$ 5%;

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






	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width: LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here

	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA

	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI Quad CRT-D SureScan, model DTMB2QQ consisting of: – implantable cardioverter defibrillator AMPLIA MRI Quad CRT-D SureScan, model DTMB2QQ – torque wrench; – Device Manual AMPLIA MRI Quad CRT-D SureScan DTMB2QQ; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliu 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA

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- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)
- tissue necrosis⁴
- valve damage (particularly in fragile hearts)

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Cardiology requirements

- The patient has no implanted lead extenders, lead adaptors, or abandoned leads.
- The patient has no broken leads or leads with intermittent electrical contact, as confirmed by lead impedance history.
- The SureScan CRT-D system is implanted in the left or right pectoral region.
- The SureScan device is operating within the projected service life.

	• For patients whose device will be programmed to an asynchronous pacing mode when MRI SureScan is programmed to On, no diaphragmatic stimulation is present when the paced leads have a pacing output of 5.0 V and a pulse width of 1.0 ms. Note: The LV lead is not paced during SureScan operation so the presence of diaphragmatic stimulation on the LV lead at a pacing output of 5.0 V and a pulse width of 1.0 ms does not need to be considered. Caution: For pacemaker-dependent patients, it is not recommended to perform an MRI scan if the right ventricular (RV) lead pacing capture threshold is greater than 2.0 V at 0.4 ms. A higher pacing capture threshold may indicate an issue with the implanted lead. Patient monitoring and rescue requirements • Continuous patient monitoring is required while MRI SureScan mode is programmed to On. • In the event that patient rescue is required, an external defibrillator must be immediately available. Note: For radiology requirements for an MRI scan, refer to the MRI technical manual. Pre-implant consideration Patient evaluation for the implant of Amplia MRI different models should include the following consideration about a concomitant implant with a neurostimulator: Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads; • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected); • The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI Quad CRT-D SureScan, model DTMB2QQ	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper

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Sterile Package:

Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm $\pm$ 1 cm x 13 cm $\pm$ 0.5 cm x 2.1 cm $\pm$ 0.4 cm;

The mass of the package with the device, g: 121 $\pm$ 5%.

Secondary Package:

Dimensions, L x W x H: 17.54 cm $\pm$ 0.8 cm x 14.31 cm $\pm$ 0.7 cm x 4.51 cm $\pm$ 0.18 cm;



The mass of the package with the device, g: 513 $\pm$ 5%.

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
Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
REF	Reorder number
SN	Serial number
	Use by
	Consult instructions for use
	Date of manufacture
	Package contents

	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here

	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA

	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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

12. Information on the process of designing, developing **TELCSECURITY_RW1_1 version 3.1 and upper** – firmware Implantable Cardioverter Defibrillator
 Release Date: 06.04.2020.
 Safety class – C (regarding IEC 62304).

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13. Sterilization method, Sterilization Validation

Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced) Routine Sterilization Parameter Specifications

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disposal or destruction of a medical device	• To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. • Explant the implanted device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. • Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. • Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. Addresses are listed below.	
	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
Waste class – B.		
15. Shelf life period	18 months.	
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.	
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.	
18. Warranty	Guaranteed shelf life -18 months;	

19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS	COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled "Sterile"- Part 1: Requirements for terminally sterilized medical devices	Full Compliance
	EN 980 Symbols for use in the labelling of medical devices	Full Compliance
	EN 1041 Information supplied by the manufacturer of medical devices	Full Compliance
	ISO 5841-3 Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
	EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
	EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
	EN ISO 11135 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
	ISO 11318 Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
	EN ISO 11607-1 Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
	EN ISO 14971 Medical devices - Application of risk management to medical devices	Full Compliance
	ISO 27186 Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
	EN 45502-1 Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
	EN 45502-2-1 Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI CRT-D SureScan, model DTMB2D4 consisting of: – implantable cardioverter defibrillator AMPLIA MRI CRT-D SureScan, model DTMB2D4 – torque wrench; – Device Manual AMPLIA MRI CRT-D SureScan DTMB2D4; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliat 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA


	Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Amplia MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device) ● air embolism ● bleeding

- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)
- tissue necrosis*
- valve damage (particularly in fragile hearts)


	• venous occlusion • venous perforation An additional potential adverse event associated with the use of transvenous left ventricular pacing leads is coronary sinus dissection. Additional known potential adverse events associated with the use of ICD systems include the following events: • inappropriate shocks • potential mortality due to inability to defibrillate • shunting current or insulating myocardium during defibrillation Patients susceptible to frequent shocks despite medical management could develop psychological intolerance to an ICD system that might include the following conditions: • dependency • depression • fear of premature battery depletion • fear of shocking while conscious • fear that shocking capability may be lost • imagined shocking (phantom shock) The SureScan CRT-D system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment: • lead electrode heating and tissue damage resulting in loss of sensing or capture or both; • spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while the MRI SureScan mode is programmed to On; • potential for ventricular tachycardia or ventricular fibrillation (VT/VF) induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan mode; • device heating resulting in tissue damage in the implant pocket or patient discomfort or both; • MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three; • damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly; • damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer; • movement or vibration of the device or leads resulting in dislodgment.
5. Condition for use	The product is intended for use in medical institutions. Use of a pacemaker is permitted only under the supervision or supervision of specially trained medical personnel. Proper surgical procedures and sterile techniques are the responsibility of the physician. Each physician must apply the information in these procedures according to professional medical training and experience. Operational conditions: internal body environment +32°C-+42°C, humidity 80-100%.
6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator. 47206, Cardiac pulse generator software.

Class depending on the potential risk of a medical product - CLASS 3 / III.
 Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs.
 Safety Class of Software in accordance with IEC 62304 – C.
 Type of contact with the body - Long-term exposure with blood, tissues and body fluids.
 Sterility: sterile.
 Classification, depending on the possible consequences of failure during the use of the product: A.

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
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

Radiopaque location and dimensions (inches(mm)):

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
- the Amplia MRI / Claria MRI / Compia MRI device
- a SureScan right atrial pacing lead or a Model 6725 pin plug for the right atrial port
- a SureScan left ventricular pacing lead
- a SureScan defibrillation lead

When a single coil SureScan defibrillation lead is used, a Medtronic DF-1 pin plug must be secured in the SVC port to make a complete SureScan defibrillation system. DF-1 pin plug only for models AMPLIA MRI Quad CRT-D SureScan DTMB2Q1, AMPLIA MRI CRT-D SureScan DTMB2D1, CLARIA MRI Quad CRT-D SureScan DTMA2Q1, CLARIA MRI CRT-D SureScan DTMA2D1, COMPIA MRI CRT-D SureScan DTMC2D1.

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Cardiology requirements

- The patient has no implanted lead extenders, lead adaptors, or abandoned leads.

	• The patient has no broken leads or leads with intermittent electrical contact, as confirmed by lead impedance history. • The SureScan CRT-D system is implanted in the left or right pectoral region. • The SureScan device is operating within the projected service life. • For patients whose device will be programmed to an asynchronous pacing mode when MRI SureScan is programmed to On, no diaphragmatic stimulation is present when the paced leads have a pacing output of 5.0 V and a pulse width of 1.0 ms. Note: The LV lead is not paced during SureScan operation so the presence of diaphragmatic stimulation on the LV lead at a pacing output of 5.0 V and a pulse width of 1.0 ms does not need to be considered. Caution: For pacemaker-dependent patients, it is not recommended to perform an MRI scan if the right ventricular (RV) lead pacing capture threshold is greater than 2.0 V at 0.4 ms. A higher pacing capture threshold may indicate an issue with the implanted lead. Patient monitoring and rescue requirements • Continuous patient monitoring is required while MRI SureScan mode is programmed to On. • In the event that patient rescue is required, an external defibrillator must be immediately available. Note: For radiology requirements for an MRI scan, refer to the MRI technical manual. Pre-implant consideration Patient evaluation for the implant of Amplia MRI / Claria MRI / Compia MRI different models should include the following consideration about a concomitant implant with a neurostimulator: Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads;

- Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected);
- The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

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Device name	Material Type
Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI CRT-D SureScan, model DTMB2D4	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V /

11. Information about Labelling and Packaging of medical device

Package content

- 1 implantable cardioverter defibrillator
- 1 torque wrench
- Device Manual;
- Limited warranty;
- Warnings and precautions related to medical procedures and EMF;
- MRI Technical Manual;
- Explanation of symbols for ICT and CRT-D devices;
- Document "Electronic operating instructions / manuals";
- CD-ROM.

Information on labelling, packaging

Sterile Package:

Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm;

The mass of the package with the device, g: 120 ± 5%;

Secondary Package:

Dimensions, L x W x H: 17.54 cm ± 0.8 cm x 14.31 cm ± 0.7 cm x 4.51 cm ± 0.18 cm;

The mass of the package with the device, g: 512 ± 5%;

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Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
	Reorder number
	Serial number
	Use by
	Consult instructions for use

	Date of manufacture
	Package contents
	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width LV
	Sensitivity: RV
	Sensitivity: RA

	A-V interval (paced and sensed)
	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion

	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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

12. Information on the process of designing, developing

TELCSECURITY_RW1_1 version 3.1 and upper – firmware Implantable Cardioverter Defibrillator
 Release Date: 06.04.2020.
 Safety class – C (regarding IEC 62304).

and validating software used in a finished medical device

13. Sterilization method, Sterilization Validation

Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)

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disposal or destruction of a medical device	• To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. • Explant the implanted device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. • Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. • Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. Addresses are listed below.	
	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
	Waste class – B.	
15. Shelf life period	18 months.	
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.	
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.	
18. Warranty	Guaranteed shelf life -18 months;	

19. The list of national and international regulatory documents/standards

REFERENCED STANDARDS & REVISIONS		COMPLIANCE LEVEL
EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled "Sterile"- Part 1: Requirements for terminally sterilized medical devices		Full Compliance
EN 980 Symbols for use in the labelling of medical devices		Full Compliance
EN 1041 Information supplied by the manufacturer of medical devices		Full Compliance
ISO 5841-3 Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers		Full Compliance
EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process		Full Compliance
EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals		Full Compliance
EN ISO 11135 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices		Full Compliance
ISO 11318 Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement		Full Compliance
EN ISO 11607-1 Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems		Full Compliance
EN ISO 14971 Medical devices - Application of risk management to medical devices		Full Compliance
ISO 27186 Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements		Full Compliance
EN 45502-1 Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer		Full Compliance
EN 45502-2-1 Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)		Full Compliance to Section 6

	EN 45502-2-2 Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
	EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
	EN 60601-1-6 Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI Quad CRT-D SureScan, model DTMC2QQ consisting of: – implantable cardioverter defibrillator COMPIA MRI Quad CRT-D SureScan, model DTMC2QQ – torque wrench; – Device Manual COMPIA MRI Quad CRT-D SureScan DTMC2QQ; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA

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
- bleeding
- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)
- tissue necrosis

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H (mm) × W (mm) × D (mm)	74± 1 mm × 51± 1 mm × 13± 1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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- The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI Quad CRT-D SureScan, model DTMC2QQ	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™

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The mass of the package with the device, g: 121 ± 5%.

Secondary Package:


Dimensions, L x W x H: 17.54 cm ± 0.8 cm x 14.31 cm ± 0.7 cm x 4.51 cm ± 0.18 cm;




The mass of the package with the device, g: 513 ± 5%.

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Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
REF	Reorder number
SN	Serial number
	Use by
	Consult instructions for use
	Date of manufacture

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

Routine Sterilization Parameter Specifications

Cycle Phase	Sterilization Parameter	Acceptance Criteria
Conditioning	Vacuum Rate	7.0 – 13.0 kPa/minute
	Pressure at the beginning of leak check	14.0 - 20.0 kPa
	Leak Check Duration	15 ± 2 minutes
	Leak Check Vacuum Loss	< 4 kPa
	Time	30 – 35 minutes
	Number of Humidity Pulses	4
	Temperature (Chamber and Load) ¹	45.0 – 55.0°C
	Chamber Relative Humidity ¹	55.0 – 85.0% RH
	Pressure, just prior to gas injection	8.0 – 14.0 kPa
	Fans (#1 and #2)	0.030 – 0.300 amps
EO Exposure	Gas Inlet Temperature	15.0 – 75.0°C
	Exposure Time	27 - 33 minutes
	Temperature (Chamber and Load) ²	45.0 – 55.0°C
	Pressure	55.0 – 75.0 kPa
	Gas Concentration ²	≥ 720 mg/L
	EO Gas Cartridge Weight	127 ± 5 grams
Post Vacuum	Fans (#1 and #2)	0.030 – 0.300 amps
	Vacuum Rate	1.0 - 2.0 kPa/minute
Air Flush	Pressure	20.0 - 30.0 kPa
	Time	≥ 25 minutes
Aeration	Temperature	44.0 – 56.0°C
	Time	MJC: 2 – 16 hours (1X sterilization); 8 – 16 hours (2X, 3X, 4X, or 5X sterilization) SMO: 2 – 16 hours (1X sterilization); 12 – 16 hours (2X, 3X, 4X, or 5X sterilization)
	Fans (#1 and #2)	0.030 – 0.300 amps

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
July 24, 2020

Date

7/24/2020**Addition to device manual***(When the device is delivered to the Russian Federation, the information contained in this addition to device manual is a priority.)*

1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI CRT-D SureScan, model DTMC2D4 consisting of: – implantable cardioverter defibrillator COMPIA MRI CRT-D SureScan, model DTMC2D4 – torque wrench; – Device Manual COMPIA MRI CRT-D SureScan DTMC2D4; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA

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
	Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Compia MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device) ● air embolism

- bleeding
- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage*due to heating of device or lead (during an MRI procedure)
- tissue necrosis

	• valve damage (particularly in fragile hearts) • venous occlusion • venous perforation An additional potential adverse event associated with the use of transvenous left ventricular pacing leads is coronary sinus dissection. Additional known potential adverse events associated with the use of ICD systems include the following events: • inappropriate shocks • potential mortality due to inability to defibrillate • shunting current or insulating myocardium during defibrillation Patients susceptible to frequent shocks despite medical management could develop psychological intolerance to an ICD system that might include the following conditions: • dependency • depression • fear of premature battery depletion • fear of shocking while conscious • fear that shocking capability may be lost • imagined shocking (phantom shock) The SureScan CRT-D system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment: • lead electrode heating and tissue damage resulting in loss of sensing or capture or both; • spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while the MRI SureScan mode is programmed to On; • potential for ventricular tachycardia or ventricular fibrillation (VT/VF) induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan mode; • device heating resulting in tissue damage in the implant pocket or patient discomfort or both; • MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three; • damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly; • damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer; • movement or vibration of the device or leads resulting in dislodgment.
5. Condition for use	The product is intended for use in medical institutions. Use of a pacemaker is permitted only under the supervision or supervision of specially trained medical personnel. Proper surgical procedures and sterile techniques are the responsibility of the physician. Each physician must apply the information in these procedures according to professional medical training and experience. Operational conditions: internal body environment +32°C-+42°C, humidity 80-100%.
6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator.

47206, Cardiac pulse generator software.
Class depending on the potential risk of a medical product - CLASS 3 / III.
Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs.
Safety Class of Software in accordance with IEC 62304 – C.
Type of contact with the body - Long-term exposure with blood, tissues and body fluids.
Sterility: sterile.
Classification, depending on the possible consequences of failure during the use of the product: A.

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
H (mm) × W (mm) × D (mm)	73± 1 mm × 51± 1 mm × 13± 1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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
The device responds to bradyarrhythmias by providing bradycardia pacing therapy. LV-only, RV-synchronized LV-only, simultaneous or sequential BiV pacing is used to provide patients with CRT. The device also provides diagnostic and

monitoring information that assists with system evaluation and patient care. The device also provides treatment with defibrillation, cardioversion, and ATP therapies. The devices can also detect atrial tachyarrhythmias (AT and AF).

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Cardiology requirements

- The patient has no implanted lead extenders, lead adaptors, or abandoned leads.
- The patient has no broken leads or leads with intermittent electrical contact, as confirmed by lead impedance history.
- The SureScan CRT-D system is implanted in the left or right pectoral region.

	• The SureScan device is operating within the projected service life. • For patients whose device will be programmed to an asynchronous pacing mode when MRI SureScan is programmed to On, no diaphragmatic stimulation is present when the paced leads have a pacing output of 5.0 V and a pulse width of 1.0 ms. Note: The LV lead is not paced during SureScan operation so the presence of diaphragmatic stimulation on the LV lead at a pacing output of 5.0 V and a pulse width of 1.0 ms does not need to be considered. Caution: For pacemaker-dependent patients, it is not recommended to perform an MRI scan if the right ventricular (RV) lead pacing capture threshold is greater than 2.0 V at 0.4 ms. A higher pacing capture threshold may indicate an issue with the implanted lead. Patient monitoring and rescue requirements • Continuous patient monitoring is required while MRI SureScan mode is programmed to On. • In the event that patient rescue is required, an external defibrillator must be immediately available. Note: For radiology requirements for an MRI scan, refer to the MRI technical manual. Pre-implant consideration Patient evaluation for the implant of Compia MRI different models should include the following consideration about a concomitant implant with a neurostimulator: Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads; • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected);

• The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI CRT-D SureScan, model DTMC2D4	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ -Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™

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Sterile Package:

Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm;
The mass of the package with the device, g: 120 ± 5%.


Secondary Package:


Dimensions, L x W x H: 17.54 cm ± 0.8 cm x 14.31 cm ± 0.7 cm x 4.51 cm ± 0.18 cm;
The mass of the package with the device, g: 512 ± 5%.

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
Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
	Reorder number
	Serial number
	Use by
	Consult instructions for use
	Date of manufacture

	Package contents
	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width: LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)

	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion

	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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disposal or destruction of a medical device	• To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. • Explant the implanted device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. • Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. • Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. Addresses are listed below.	
	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliat 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
Waste class – B.		
15. Shelf life period	18 months.	
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.	
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.	
18. Warranty	Guaranteed shelf life -18 months;	

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	Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
	EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
	EN 60601-1-6 Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI Quad CRT-D SureScan, model DTMB2Q1 consisting of: – implantable cardioverter defibrillator AMPLIA MRI Quad CRT-D SureScan, model DTMB2Q1; – DF-1 pin plug; – torque wrench; – Device Manual AMPLIA MRI Quad CRT-D SureScan DTMB2Q1; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

	3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Amplia MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device)

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6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator. 47206, Cardiac pulse generator software. Class depending on the potential risk of a medical product - CLASS 3 / III. Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs. Safety Class of Software in accordance with IEC 62304 – C. Type of contact with the body - Long-term exposure with blood, tissues and body fluids. Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A.
7. Design and appearance, photo or image, the General Scheme (flowchart)	AMPLIA MRI Quad CRT-D DTMB2Q1 Torque wrench DF-1 Connector Pin Plug
8. Description of the main functional parts of a medical	Physical characteristics of implantable device:

device and technical details.

Specification	DTMB2Q1
Volume	36 cm ³ ±0.5 cm ³
Mass	82 g± 1 g
H (mm) × W (mm) × D (mm)	74±1 mm × 51±1 mm × 13±1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.

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
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
11. Information about Labelling and Packaging of medical device	Package content - 1 implantable cardioverter defibrillator - 1 torque wrench - 1 DF-1 pin plug - Device Manual; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM. Information on 'labelling, packaging Sterile Package:	







Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm;
The mass of the package with the device, g: 122 ± 5%.
Secondary Package:
Dimensions, L x W x H: 17.54 cm ± 0.8 cm x 14.31 cm ± 0.7 cm x 4.51 cm ± 0.18 cm;
The mass of the package with the device, g: 514 ± 5%.

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	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width: LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT

	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional

	Symbol SureScan
DDD	DDD

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

	
12. Information on the process of designing, developing and validating software used in a finished medical device	TELCSECURITY_RW1_1 version 3.1 and upper – firmware Implantable Cardioverter Defibrillator Release Date: 06.04.2020. Safety class – C (regarding IEC 62304).
13. Sterilization method, Sterilization Validation	Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced) Routine Sterilization Parameter Specifications

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	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com	
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000	
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800	
	Waste class – B.		
15. Shelf life period	18 months.		
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.		
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.		
18. Warranty	Guaranteed shelf life -18 months; Warranty period - 5 years.		
19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS		COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled “Sterile”- Part 1: Requirements for terminally sterilized medical devices		Full Compliance
	EN 980 Symbols for use in the labelling of medical devices		Full Compliance

EN 1041	Information supplied by the manufacturer of medical devices	Full Compliance
ISO 5841-3	Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
EN ISO 10993-7	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
EN ISO 11135	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
ISO 11318	Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
EN ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
EN ISO 14971	Medical devices - Application of risk management to medical devices	Full Compliance
ISO 27186	Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
EN 45502-1	Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
EN 45502-2-1	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6
EN 45502-2-2	Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
EN 60601-1	Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
EN 60601-1-6 *		Full Compliance

	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1	Full Compliance
	Medical device-Application of usability engineering to medical devices	
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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

Minori Carlsson
Senior Regulatory Affairs Director
Medtronic, Inc.

July 24, 2020
Date

710 Medtronic Parkway
Minneapolis, MN 55432
USA

Cindi A Bigelbach
7/24/2020

Addition to device manual

(When the device is delivered to the Russian Federation, the information contained in this addition to device manual is a priority.)

1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI CRT-D SureScan, model DTMB2D1 consisting of: – implantable cardioverter defibrillator AMPLIA MRI CRT-D SureScan, model DTMB2D1 – DF-1 pin plug; – torque wrench; – Device Manual AMPLIA MRI CRT-D SureScan DTMB2D1; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

	3. Medtronic, Inc., 8200 Coral Sea Street N.E. Nounds View, MN 55112, USA Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Amplia MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: • Patients with ventricular dyssynchrony • Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block • Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: • concomitant implant of a unipolar pacemaker • incessant VT or VF • a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF • tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: • acceleration of tachyarrhythmias (caused by device)

- air embolism
- bleeding
- body rejection phenomena, including local tissue reaction
- cardiac dissection
- cardiac perforation
- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)

- tissue necrosis
- valve damage (particularly in fragile hearts)
- venous occlusion
- venous perforation

An additional potential adverse event associated with the use of transvenous left ventricular pacing leads is coronary sinus dissection.

Additional known potential adverse events associated with the use of ICD systems include the following events:

- inappropriate shocks
- potential mortality due to inability to defibrillate
- shunting current or insulating myocardium during defibrillation

Patients susceptible to frequent shocks despite medical management could develop psychological intolerance to an ICD system that might include the following conditions:

- dependency
- depression
- fear of premature battery depletion
- fear of shocking while conscious
- fear that shocking capability may be lost
- imagined shocking (phantom shock)

The SureScan CRT-D system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment:

- lead electrode heating and tissue damage resulting in loss of sensing or capture or both;
- spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while the MRI SureScan mode is programmed to On;
- potential for ventricular tachycardia or ventricular fibrillation (VT/VF) induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan mode;
- device heating resulting in tissue damage in the implant pocket or patient discomfort or both;
- MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three;
- damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly;
- damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer;
- movement or vibration of the device or leads resulting in dislodgment.

5. Condition for use

The product is intended for use in medical institutions. Use of a pacemaker is permitted only under the supervision or supervision of specially trained medical personnel.


Proper surgical procedures and sterile techniques are the responsibility of the physician. Each physician must apply the information in these procedures according to professional medical training and experience. Operational conditions: internal body environment +32°C-+42°C, humidity 80-100%.

6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator. 47206, Cardiac pulse generator software. Class depending on the potential risk of a medical product - CLASS 3 / III. Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs. Safety Class of Software in accordance with IEC 62304 – C. Type of contact with the body - Long-term exposure with blood, tissues and body fluids. Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A.
7. Design and appearance, photo or image, the General Scheme (flowchart)	AMPLIA MRI CRT-D DTMB2D1 Torque wrench DF-1 Connector Pin Plug
8. Description of the main functional	Physical characteristics of implantable device:

parts of a medical device and technical details.

Specification	DTMB2D1
Volume	35 cm ³ ±0.5 cm ³
Mass	80 g± 1 g
H (mm) × W (mm) × D (mm)	71± 1 mm x 51± 1 mm x 13± 1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

Radiopaq location and dimensions (inches(mm)):

0,008 in = 0,203 mm, 0,090 in = 2,286 mm, 0,310 in = 7,874 mm

Principles of operation

The basic functions of the CRT-D devices are to provide BiV or LV pacing for CRT, to monitor and regulate a patient's heart rate by providing dual chamber rate-responsive bradycardia pacing, to provide diagnostics, and to provide treatment for atrial and ventricular arrhythmias. The lead system used with these devices must provide pacing to the LV, sensing, pacing, and cardioversion/defibrillation therapies to the RV, and sensing and pacing to the atrium. The CRTD device and the pacing leads constitute the implantable portion of the device system. The device senses the electrical activity of the patient's heart using the electrodes of the implanted atrial and RV leads. It then analyzes the heart rhythm based on selectable


detection parameters. The device monitors the heart rhythm for ventricular tachyarrhythmias (VT / VF) and uses detection criteria to distinguish between true ventricular arrhythmias and rapidly conducted SVT.

The device responds to bradyarrhythmias by providing bradycardia pacing therapy. LV-only, RV-synchronized LV-only, simultaneous or sequential BiV pacing is used to provide patients with CRT. The device also provides diagnostic and monitoring information that assists with system evaluation and patient care. The device also provides treatment with defibrillation, cardioversion, and ATP therapies. The devices can also detect atrial tachyarrhythmias (AT and AF).

Technical details for accessories supplied with main device:

Torque wrench:

Torque wrench – Use only the torque wrench supplied with the device. The torque wrench is designed to prevent damage to the device from overtightening a setscrew. Other torque wrenches (for example, a blue-handled or right-angled hex wrench) have torque capabilities greater than the lead connector can tolerate.


The torque is limited to 12 to 17 in-oz. during first clockwise revolution; it will click when torque value in this range is achieved

Mass (g): 8,0±0,5.

DF-1 connector pin plug:

DF-1 Connector Pin Plug (unit: mm)

*

This device is a plug to fill an IS-1 (ISO 5841-3) compliant connector port of pulse generators, such as an implantable cardiac pacemaker, to which an implantable IPG/ICD lead connector is not connected.

Mass (g): 0.41 ± 0.005.

MRI conditions for use

A complete SureScan CRT-D system is required for use in the MR environment. A complete SureScan CRT-D system includes the following components. Any other combination may result in a hazard to the patient during an MRI scan. To verify that components are part of a SureScan system, visit <http://www.mrisurescan.com>.

- the Amplia MRI device
- a SureScan right atrial pacing lead or a Model 6725 pin plug for the right atrial port
- a SureScan left ventricular pacing lead
- a SureScan defibrillation lead

When a single coil SureScan defibrillation lead is used, a Medtronic DF-1 pin plug must be secured in the SVC port to make a complete SureScan defibrillation system.

Warning: Do not scan a patient without first programming MRI SureScan to On. Scanning the patient without programming MRI SureScan to On may result in patient harm or damage to the SureScan CRT-D system.

Note: MRI SureScan cannot be programmed to On if the device is recommended for replacement.

Patients and their implanted systems must be screened to meet the following requirements:

Cardiology requirements

- The patient has no implanted lead extenders, lead adaptors, or abandoned leads.
- The patient has no broken leads or leads with intermittent electrical contact, as confirmed by lead impedance history.
- The SureScan CRT-D system is implanted in the left or right pectoral region.
- The SureScan device is operating within the projected service life.
- For patients whose device will be programmed to an asynchronous pacing mode when MRI SureScan is programmed to On, no diaphragmatic stimulation is present when the paced leads have a pacing output of 5.0 V and a pulse width of 1.0 ms.

Note: The LV lead is not paced during SureScan operation so the presence of diaphragmatic stimulation on the LV lead at a pacing output of 5.0 V and a pulse width of 1.0 ms does not need to be considered.

Caution: For pacemaker-dependent patients, it is not recommended to perform an MRI scan if the right ventricular (RV) lead pacing capture threshold is greater than 2.0 V at 0.4 ms. A higher pacing capture threshold may indicate an issue with the implanted lead.

Patient monitoring and rescue requirements

- Continuous patient monitoring is required while MRI SureScan mode is programmed to On.
- In the event that patient rescue is required, an external defibrillator must be immediately available.

Note: For radiology requirements for an MRI scan, refer to the MRI technical manual.

Pre-implant consideration

Patient evaluation for the implant of Amplia MRI different models should include the following consideration about a concomitant implant with a neurostimulator:

Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices

	that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.	
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads; • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected); • The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.	
10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology AMPLIA MRI CRT-D SureScan, model DTMB2D1	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Pin plug Medtronic DF-1:	Handle - liquid silicone rubber (HSC) of the biomedical class "Silastic" Stem - 316 L stainless steel

	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
11. Information about Labelling and Packaging of medical device Package content - 1 implantable cardioverter defibrillator - 1 torque wrench - 1 DF-1 pin plug - Device Manual; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM. Information on labelling, packaging Sterile Package:		

Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm $\pm$ 1 cm x 13 cm $\pm$ 0.5 cm x 2.1 cm $\pm$ 0.4 cm;
The mass of the package with the device, g: 120 $\pm$ 5%.



Secondary Package:

Dimensions, L x W x H: 17.54 cm $\pm$ 0.8 cm x 14.31 cm $\pm$ 0.7 cm x 4.51 cm $\pm$ 0.18 cm;





The mass of the package with the device, g: 512 $\pm$ 5%.

Explanation of symbols

Refer to the product package labels to see which symbols apply to this product.

Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
	Reorder number
	Serial number
	Use by
	Consult instructions for use
	Date of manufacture
	Package contents
	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage

	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT

	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional

	Symbol SureScan
DDD	DDD


QR-code


ME PRODUCTS AND ME SYSTEMS, including radio frequency transmitters or intentionally using radio frequency electromagnetic energy for diagnosis or treatment, should be marked as a sign of non-ionizing radiation

Example of device labeling (front side) – labeling can be different depend on model number of device:

Example of device labeling (back side) - labeling can be different depend on model number of device:

	 Connector port LV RV ACTIVE CAN ON/OFF Barcode IS-1 DF4-LLHH Description of connector type
12. Information on the process of designing, developing and validating software used in a finished medical device	TELCSECURITY_RW1_1 version 3.1 and upper – firmware Implantable Cardioverter Defibrillator Release Date: 06.04.2020. Safety class – C (regarding IEC 62304).
13. Sterilization method, Sterilization Validation	Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced) Routine Sterilization Parameter Specifications

Cycle Phase	Sterilization Parameter	Acceptance Criteria
Conditioning	Vacuum Rate	7.0 – 13.0 kPa/minute
	Pressure at the beginning of leak check	14.0 - 20.0 kPa
	Leak Check Duration	15 ± 2 minutes
	Leak Check Vacuum Loss	< 4 kPa
	Time	30 – 35 minutes
	Number of Humidity Pulses	4
	Temperature (Chamber and Load) ¹	45.0 – 55.0°C
	Chamber Relative Humidity ¹	55.0 – 85.0% RH
	Pressure, just prior to gas injection	8.0 – 14.0 kPa
	Fans (#1 and #2)	0.030 – 0.300 amps
EO Exposure	Gas Inlet Temperature	15.0 – 75.0°C
	Exposure Time	27 - 33 minutes
	Temperature (Chamber and Load) ²	45.0 – 55.0°C
	Pressure	55.0 – 75.0 kPa
	Gas Concentration ²	≥ 720 mg/L
	EO Gas Cartridge Weight	127 ± 5 grams
	Fans (#1 and #2)	0.030 – 0.300 amps
Post Vacuum	Vacuum Rate	1.0 - 2.0 kPa/minute
	Pressure	20.0 - 30.0 kPa
Air Flush	Time	≥ 25 minutes
Aeration	Temperature	44.0 – 56.0°C
	Time	MJC: 2 – 16 hours (1X sterilization); 8 – 16 hours (2X, 3X, 4X, or 5X sterilization) SMO: 2 – 16 hours (1X sterilization); 12 – 16 hours (2X, 3X, 4X, or 5X sterilization)
	Fans (#1 and #2)	0.030 – 0.300 amps

¹ Relative humidity & temperature requirements must be met at the end of conditioning – just prior to gas inject.

² During the first five minutes of exposure (directly after the completion of gas injection), temperature and gas concentration spikes or drops may exceed the process parameters while the chamber is equilibrating.

All internal and external requirements have been met. The CRT-D MRI product family is qualified for 1X – 5X sterilization in the 30-minute 100% EO sterilization process. Sterility assurance level (SAL) of 10⁻⁶.

14. Procedure and conditions for the disposal or destruction of a medical device

Consider the following information related to device explant and disposal:

- To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device.
- Explant the implanted device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures.
- Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices.
- Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. Addresses are listed below.

	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
	Waste class – B.	
15. Shelf life period	18 months.	
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.	
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.	
18. Warranty	Guaranteed shelf life -18 months; Warranty period - 5 years.	
19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS	COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled “Sterile”- Part 1: Requirements for terminally sterilized medical devices	Full Compliance
	EN 980 Symbols for use in the labelling of medical devices	Full Compliance

	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
	Waste class – B.	
15. Shelf life period	18 months.	
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.	
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.	
18. Warranty	Guaranteed shelf life -18 months; Warranty period - 5 years.	
19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS	COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled “Sterile”- Part 1: Requirements for terminally sterilized medical devices	Full Compliance
	EN 980 Symbols for use in the labelling of medical devices	Full Compliance

EN 1041 Information supplied by the manufacturer of medical devices	Full Compliance
ISO 5841-3 Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
EN ISO 11135 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
ISO 11318 Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
EN ISO 11607-1 Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
EN ISO 14971 Medical devices - Application of risk management to medical devices	Full Compliance
ISO 27186 Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
EN 45502-1 Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
EN 45502-2-1 Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6
EN 45502-2-2 Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
EN 60601-1-6*	Full Compliance

	Medical electrical equipment -- Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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
[Handwritten signature]

Minori Carlsson
Minori Carlsson
Senior Regulatory Affairs Director
Medtronic, Inc.

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI Quad CRT-D SureScan, model DTMA2Q1 consisting of: – implantable cardioverter defibrillator CLARIA MRI Quad CRT-D SureScan, model DTMA2Q1 – DF-1 pin plug; – torque wrench; – Device Manual CLARIA MRI Quad CRT-D SureScan DTMA2Q1; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

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	3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Claria MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device)

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
	● tissue necrosis ● valve damage (particularly in fragile hearts) ● venous occlusion ● venous perforation An additional potential adverse event associated with the use of transvenous left ventricular pacing leads is coronary sinus dissection. Additional known potential adverse events associated with the use of ICD systems include the following events: ● inappropriate shocks ● potential mortality due to inability to defibrillate ● shunting current or insulating myocardium during defibrillation Patients susceptible to frequent shocks despite medical management could develop psychological intolerance to an ICD system that might include the following conditions: ● dependency ● depression ● fear of premature battery depletion ● fear of shocking while conscious ● fear that shocking capability may be lost ● imagined shocking (phantom shock) The SureScan CRT-D system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment: ● lead electrode heating and tissue damage resulting in loss of sensing or capture or both; ● spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while the MRI SureScan mode is programmed to On; ● potential for ventricular tachycardia or ventricular fibrillation (VT/VF) induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan mode; ● device heating resulting in tissue damage in the implant pocket or patient discomfort or both; ● MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three; ● damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly; ● damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer; ● movement or vibration of the device or leads resulting in dislodgment.
5. Condition for use	The product is intended for use in medical institutions. Use of a pacemaker is permitted only under the supervision or supervision of specially trained medical personnel. Proper surgical procedures and sterile techniques are the responsibility of the physician. Each physician must apply the information in these procedures according to professional medical training and experience. Operational conditions: internal body environment +32°C-+42°C, humidity 80-100%.


6. Classification of medical product:

GMDN
47270, Cardiac resynchronization therapy implantable defibrillator.
47206, Cardiac pulse generator software.
Class depending on the potential risk of a medical product - CLASS 3 / III.
Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs.
Safety Class of Software in accordance with IEC 62304 – C.
Type of contact with the body - Long-term exposure with blood, tissues and body fluids.
Sterility: sterile.
Classification, depending on the possible consequences of failure during the use of the product: A.


7. Design and appearance, photo or image, the General Scheme (flowchart)

CLARIA MRI Quad CRT-D DTMA2Q1

Torque wrench

DF-1 Connector Pin Plug


8. Description of the main functional

Physical characteristics of implantable device:	
Specification	DTMA2Q1

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Volume	36 cm ³ ±0.5 cm ³
Mass	82 g± 1 g
H (mm) × W (mm) × D (mm)	74±1 mm × 51±1 mm × 13±1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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0,008 in = 0,203 mm, 0,090 in = 2,286 mm, 0,310 in = 7,874 mm

Principles of operation


The basic functions of the CRT-D devices are to provide BiV or LV pacing for CRT, to monitor and regulate a patient's heart rate by providing dual chamber rate-responsive bradycardia pacing, to provide diagnostics, and to provide treatment for atrial and ventricular arrhythmias. The lead system used with these devices must provide pacing to the LV, sensing, pacing, and cardioversion/defibrillation therapies to the RV, and sensing and pacing to the atrium. The CRTD device and the pacing leads constitute the implantable portion of the device system. The device senses the electrical activity of the patient's heart using the electrodes of the implanted atrial and RV leads. It then analyzes the heart rhythm based on selectable detection parameters. The device monitors the heart rhythm for ventricular tachyarrhythmias (VT / VF) and uses detection criteria to distinguish between true ventricular arrhythmias and rapidly conducted SVT.

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
- the Claria MRI device
- a SureScan right atrial pacing lead or a Model 6725 pin plug for the right atrial port
- a SureScan left ventricular pacing lead
- a SureScan defibrillation lead

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Pre-implant consideration

Patient evaluation for the implant of Claria MRI different models should include the following consideration about a concomitant implant with a neurostimulator:

Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant

	procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.	
9. The possibility and ways of integrating with other medical device	The CRT-D MRI System may be comprised of the following components: • CRT-D bipolar and quadripolar devices; • SW034 Programmer Software Application; • CareLink Programmer Model 2090; • CareLink Encore Programmer Model 29901; • Conexus Activator Model 27901; • CareLink Monitor Model 2490C; • MyCareLink Monitor Model 24950; • 2020B CareLink Express Monitor; • Commercially available MRI SureScan Right Atrial Pacing Leads; • Commercially available MRI SureScan Right Ventricular Defibrillation Leads; • Commercially available MRI SureScan Left Ventricular Pacing Leads; • Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected); • The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.	
10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)	Device name	Material Type
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI Quad CRT-D SureScan, model DTMA2Q1	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
	Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
	Pin plug Medtronic DF-1:	Handle - liquid silicone rubber (HSC) of the biomedical class "Silastic" Stem - 316 L stainless steel
	Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated

		Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
11. Information about Labelling and Packaging of medical device	Package content - 1 implantable cardioverter defibrillator - 1 torque wrench - 1 DF-1 pin plug - Device Manual; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM. Information on labelling, packaging Sterile Package: Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm; The mass of the package with the device, g: 122 ± 5%.	

Secondary Package:

Dimensions, L x W x H: 17.54 cm $\pm$ 0.8 cm x 14.31 cm $\pm$ 0.7 cm x 4.51 cm $\pm$ 0.18 cm;

The mass of the package with the device, g: 514 $\pm$ 5%.

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Symbols	Designation
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
REF	Reorder number
SN	Serial number
	Use by
	Consult instructions for use
	Date of manufacture
	Package contents
	Manufacturer
EC REP	Authorized representative in the European community
STERILE EO	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape

	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate
	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width: LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can

	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation
	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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

12. Information on the process of designing, developing and validating software used in a finished medical device	TELCSECURITY_RW1_1 version 3.1 and upper – firmware Implantable Cardioverter Defibrillator Release Date: 06.04.2020. Safety class – C (regarding IEC 62304).
13. Sterilization method, Sterilization Validation	Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced) Routine Sterilization Parameter Specifications

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EN 1041	Information supplied by the manufacturer of medical devices	Full Compliance
ISO 5841-3	Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
EN ISO 10993-1	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
EN ISO 10993-7	Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
EN ISO 11135	Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
ISO 11318	Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
EN ISO 11607-1	Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
EN ISO 14971	Medical devices - Application of risk management to medical devices	Full Compliance
ISO 27186	Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
EN 45502-1	Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
EN 45502-2-1	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6
EN 45502-2-2	Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
EN 60601-1	Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
EN 60601-1-6		Full Compliance

	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability			
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance		
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance		
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. <table><tr><td>Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855</td><td>Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com</td></tr></table>		Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855	Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu , www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com
Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855	Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu , www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com			

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CINDI A. BIGELBACH
NOTARY PUBLIC
MINNESOTA
My Commission Expires January 31, 2022

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology
	Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI CRT-D SureScan, model DTMA2D1 consisting of: - implantable cardioverter defibrillator CLARIA MRI CRT-D SureScan, model DTMA2D1 - DF-1 pin plug; - torque wrench; - Device Manual CLARIA MRI CRT-D SureScan DTMA2D1; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliau 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA 3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC

	123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The device is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device) ● air embolism ● bleeding ● body rejection phenomena, including local tissue reaction ● cardiac dissection ● cardiac perforation

- cardiac tamponade
- chronic nerve damage
- constrictive pericarditis
- death
- device migration
- endocarditis
- erosion
- excessive fibrotic tissue growth
- extrusion
- fibrillation or other arrhythmias
- fluid accumulation
- formation of hematomas/seromas or cysts
- heart block
- heart wall or vein wall rupture
- hemothorax
- infection
- keloid formation
- lead abrasion and discontinuity
- lead migration/dislodgment
- complications and mortality due to inability to deliver appropriate and intended therapy
- muscle and/or nerve stimulation
- myocardial damage
- myocardial irritability
- myopotential sensing
- pericardial effusion
- pericardial rub
- pneumothorax
- poor connection of the lead to the device, which may lead to oversensing, undersensing, or a loss of therapy
- stroke
- threshold elevation
- thrombotic embolism
- thrombosis
- tissue damage due to heating of device or lead (during an MRI procedure)
- tissue necrosis
- valve damage (particularly in fragile hearts)
- venous occlusion
- venous perforation

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
	Additional known potential adverse events associated with the use of ICD systems include the following events: • inappropriate shocks • potential mortality due to inability to defibrillate • shunting current or insulating myocardium during defibrillation Patients susceptible to frequent shocks despite medical management could develop psychological intolerance to an ICD system that might include the following conditions: • dependency • depression • fear of premature battery depletion • fear of shocking while conscious • fear that shocking capability may be lost • imagined shocking (phantom shock) The SureScan CRT-D system has been designed to minimize the potential adverse events that may cause patient harm. The following potential adverse events may occur in the MR environment: • lead electrode heating and tissue damage resulting in loss of sensing or capture or both • spontaneous tachyarrhythmia occurring during the scan that is not detected and treated because tachyarrhythmia detection is suspended while the MRI SureScan mode is programmed to On • potential for ventricular tachycardia or ventricular fibrillation (VT/VF) induction when the patient is programmed to an asynchronous pacing mode during MRI SureScan mode • device heating resulting in tissue damage in the implant pocket or patient discomfort or both • MR-induced stimulation on leads resulting in continuous capture, VT/VF, hemodynamic collapse, or all three • damage to the device or leads causing the system to fail to detect or treat irregular heartbeats, or causing the system to treat the patient's condition incorrectly • damage to the functionality or mechanical integrity of the device resulting in the inability of the device to communicate with the programmer • movement or vibration of the device or leads resulting in dislodgment
5. Condition for use	The product is intended for use in medical institutions. Use of a pacemaker is permitted only under the supervision or supervision of specially trained medical personnel. Proper surgical procedures and sterile techniques are the responsibility of the physician. Each physician must apply the information in these procedures according to professional medical training and experience. Operational conditions: internal body environment +32°C-+42°C, humidity 80-100%.
6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator. 47206, Cardiac pulse generator software. Class depending on the potential risk of a medical product - CLASS 3 / III. Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs. Safety Class of Software in accordance with IEC 62304 – C. Type of contact with the body - Long-term exposure contact with blood, tissues and body fluids.


Sterility: sterile.

Classification, depending on the possible consequences of failure during the use of the product: A.

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
Torque wrench

DF-1 Connector Pin Plug


8. Description of the main functional parts of a medical device and technical details.

Physical characteristics of implantable device:

Specification	DTMA2D1
Volume	35 cm ³ ±0.5 cm ³
Mass	80 g± 1 g
H (mm) × W (mm) × D (mm)	71± 1 mm x 51± 1 mm x 13± 1 mm

Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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

The torque is limited to 12 to 17 in-oz. during first clockwise revolution; it will click when torque value in this range is achieved.

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When a single coil SureScan defibrillation lead is used, a Medtronic DF-T pin plug must be secured in the SVC port to make a complete SureScan defibrillation system.

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Concomitant neurostimulator and cardiac device implants – Some patients have medical conditions that require the implant of both a neurostimulator and a cardiac device (for example, pacemaker, defibrillator, or monitor). In this case, physicians (for example, neurologist, neurosurgeon, cardiologist, and cardiac surgeon) involved with either device should contact their Medtronic representative before implanting the patient with the second device. Based on the particular devices that the physicians have prescribed, Medtronic can provide the necessary precautions and warnings related to the implant procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.

9. The possibility and ways of integrating with other medical device

The CRT-D MRI System may be comprised of the following components:


- CRT-D bipolar and quadripolar devices;
- SW034 Programmer Software Application;
- CareLink™ Programmer Model 2090;
- CareLink Encore™ Programmer Model 29901;

- Conexus Activator Model 27901;
- CareLink Monitor Model 2490C;
- MyCareLink Monitor Model 24950;
- 2020B CareLink Express Monitor;
- Commercially available MRI SureScan Right Atrial Pacing Leads;
- Commercially available MRI SureScan Right Ventricular Defibrillation Leads;
- Commercially available MRI SureScan Left Ventricular Pacing Leads;
- Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected);
- The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

10. List and description of materials of the medical device entering into direct or indirect contact with the patient's body (human body)

Device name	Material Type
Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology CLARIA MRI CRT-D SureScan, model DTMA2D1	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
Torque wrench	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
Pin plug Medtronic DF-1	Handle - liquid silicone rubber (HSC) of the biomedical class "Silastic" Stem - 316 L stainless steel
Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive

	Secondary sterile packaging Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive				
11. Information about Labelling and Packaging of medical device	Package content - 1 implantable cardioverter defibrillator - 1 torque wrench - 1 DF-1 pin plug - Device Manual; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM. Information on labelling, packaging Sterile Package: Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm; The mass of the package with the device, g: 120 ± 5%; Secondary Package: Dimensions, L x W x H: 17.54 cm ± 0.8 cm x 14.31 cm ± 0.7 cm x 4.51 cm ± 0.18 cm; The mass of the package with the device, g: 512 ± 5%. Explanation of symbols Refer to the product package labels to see which symbols apply to this product. <table> <tr> <th>Symbols</th><th>Designation</th></tr> <tr> <td>CE 0123</td><td>Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.</td></tr> </table>	Symbols	Designation	CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.
Symbols	Designation				
CE 0123	Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC.				

	This symbol means that the device fully complies with the Australian (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management Standards for radio communications products
	Reorder number
	Serial number
	Use by
	Consult instructions for use
	Date of manufacture
	Package contents
	Manufacturer
	Authorized representative in the European community
	Sterilized using ethylene oxide
	Temperature limitation
	Do not reuse
	Do not use if package is damaged
	Product documentation
	Dangerous voltage
	Cardiac resynchronization therapy with defibrillation (CRT-D)
	Header shape
	DF-1 connector
	DF4 connector
	IS4 connector
	IS-1 connector
	Adaptive
	Upper tracking rate and lower rate

	Amplitude and pulse width: RA
	Amplitude and pulse width: RV
	Amplitude and pulse width LV
	Sensitivity: RV
	Sensitivity: RA
	A-V interval (paced and sensed)
	Open here
	Implantable device
	Torque wrench
	Pin plug
	V pacing and V-V pace delay
	Detection: VT, FVT, VF
	Therapies (RV): VT, FVT
	Therapies (delivered and stored): VF
	Active Can
	Pace polarity: LV
	Therapies: AT, AF
	Intervention: AT, AF
	Defibrillation

	Burst: RV, LV
	AV cardioversion
	V cardioversion
	Burst+: RA
	Ramp: RA, RV, LV
	Ramp+: RV, LV
	Ramp: RV, LV
	50 Hz Burst RA
	MR Conditional
	Symbol SureScan
DDD	DDD

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

12. Information on the process of designing, developing and validating software used in a finished medical device

TELCSECURITY_RW1_1 version 3.1 and upper – firmware Implantable Cardioverter Defibrillator
 Release Date: 06.04.2020.
 Safety class – C (regarding IEC 62304).

13. Sterilization method, Sterilization Validation	Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)		
	Routine Sterilization Parameter Specifications		
	Cycle Phase	Sterilization Parameter	Acceptance Criteria
	Conditioning	Vacuum Rate	7.0 – 13.0 kPa/minute
		Pressure at the beginning of leak check	14.0 - 20.0 kPa
		Leak Check Duration	15 ± 2 minutes
		Leak Check Vacuum Loss	< 4 kPa
		Time	30 – 35 minutes
		Number of Humidity Pulses	4
		Temperature (Chamber and Load) ¹	45.0 – 55.0°C
		Chamber Relative Humidity ¹	55.0 – 85.0% RH
		Pressure, just prior to gas injection	8.0 – 14.0 kPa
		Fans (#1 and #2)	0.030 – 0.300 amps
	EO Exposure	Gas Inlet Temperature	15.0 – 75.0°C
		Exposure Time	27 - 33 minutes
		Temperature (Chamber and Load) ²	45.0 – 55.0°C
		Pressure	55.0 – 75.0 kPa
		Gas Concentration ²	≥ 720 mg/L
		EO Gas Cartridge Weight	127 ± 5 grams
		Fans (#1 and #2)	0.030 – 0.300 amps
	Post Vacuum	Vacuum Rate	1.0 - 2.0 kPa/minute
		Pressure	20.0 - 30.0 kPa
	Air Flush	Time	≥ 25 minutes
Aeration	Temperature	44.0 – 56.0°C	
	Time	MJC: 2 – 16 hours (1X sterilization); 8 – 16 hours (2X, 3X, 4X, or 5X sterilization) SMO: 2 – 16 hours (1X sterilization); 12 – 16 hours (2X, 3X, 4X, or 5X sterilization)	
	Fans (#1 and #2)	0.030 – 0.300 amps	
¹ Relative humidity & temperature requirements must be met at the end of conditioning – just prior to gas inject.			
² During the first five minutes of exposure (directly after the completion of gas injection), temperature and gas concentration spikes or drops may exceed the process parameters while the chamber is equilibrating.			
All internal and external requirements have been met. The CRT-D MRI product family is qualified for 1X – 5X sterilization in the 30-minute 100% EO sterilization process. Sterility assurance level (SAL) of 10 ⁻⁶ .			
14. Procedure and conditions for the disposal or destruction of a medical device	Consider the following information related to device explant and disposal:		
	● To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. ● Explant the implanted device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. ● Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices.		

• Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. Addresses are listed below.

Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000	Representative office in Moscow: LLC Medtronic Russian Federation, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000	Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia	Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800

Waste class – B.

15. Shelf life period	18 months.		
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.		
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.		
18. Warranty	Guaranteed shelf life -18 months; Warranty period - 5 years.		
19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS		COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled “Sterile”- Part 1: Requirements for terminally sterilized medical devices		Full Compliance

EN 980 Symbols for use in the labelling of medical devices	Full Compliance
EN 1041 Information supplied by the manufacturer of medical devices	Full Compliance
ISO 5841-3 Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
EN ISO 11135 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
ISO 11318 Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
EN ISO 11607-1 Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
EN ISO 14971 Medical devices - Application of risk management to medical devices	Full Compliance
ISO 27186 Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
EN 45502-1 Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
EN 45502-2-1 Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6
EN 45502-2-2 Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14

	EN 60601-1-6 Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	Full Compliance
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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1. Name of medical product	Digital implantable cardioverter defibrillator Amplia, Claria, Compia with cardiac resynchronization therapy and SureScan Technology Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI CRT-D SureScan, model DTMC2D1 consisting of: – implantable cardioverter defibrillator COMPIA MRI CRT-D SureScan, model DTMC2D1 – DF-1 pin plug; – torque wrench; – Device Manual COMPIA MRI CRT-D SureScan DTMC2D1; – Limited warranty; – Warnings and precautions related to medical procedures and EMF; – MRI Technical Manual; – Explanation of symbols for ICT and CRT-D devices; – Document "Electronic operating instructions / manuals"; – CD-ROM.
2. Manufacturer, Designer, Manufacturing sites, Authorized Representative of Manufacturer	Manufacturer: Medtronic Inc., 710 Medtronic Parkway N.E. Minneapolis, MN 55432, USA Design facility: Medtronic Inc., 8200 Coral Sea Street, Mounds View, MN 55112, USA www.medtronic.com Addresses of manufacturing facility of medical device 1. Medtronic Europe Sarl, Route Du Molliat 31, Case Postale, 1131 Tolochenaz, SWITZERLAND 2. Medtronic Puerto Rico Operations Co. Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos PR 00777 USA

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	3. Medtronic, Inc., 8200 Coral Sea Street N.E. Mounds View, MN 55112, USA Authorized Representative of Manufacturer (Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, -mail, website): Medtronic LLC 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Tel.: (495) 580-73-77, Fax: (495) 580-73-78, www.medtronic.ru
3. Description of the medical device	The device is intended to provide atrial and/or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial and/or life-threatening ventricular tachyarrhythmias. Application: Cardiology.
4. Indications, Contraindications and Adverse events	Indications of Use: The system Compia MRI is indicated for use in patients who are at high risk of sudden death due to ventricular tachyarrhythmias and who have any of the following types of heart failure: ● Patients with ventricular dyssynchrony ● Patients with reduced ejection fraction (EF), regardless of New York Heart Association (NYHA) class, who are indicated for ventricular pacing, and have high degree AV block ● Patients with reduced EF and all of the following characteristics: – Have received a conventional pacemaker or an ICD – Have subsequently developed worsening HF despite optimal medical therapy – Have a high proportion of RV pacing The device is intended to provide atrial or ventricular antitachycardia pacing, cardioversion, and defibrillation for automated treatment of atrial or life-threatening ventricular tachyarrhythmias. Contraindications: Contraindicated for patients with the following: ● concomitant implant of a unipolar pacemaker ● incessant VT or VF ● a primary disorder of chronic atrial tachyarrhythmia with no concomitant VT or VF ● tachyarrhythmias with transient or reversible causes, including the following known issues: – acute myocardial infarction – drug intoxication – drowning – electric shock – electrolyte imbalance – hypoxia – sepsis Adverse events The following are known potential adverse events associated with the use of pacing systems: ● acceleration of tachyarrhythmias (caused by device)

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



6. Classification of medical product:	GMDN 47270, Cardiac resynchronization therapy implantable defibrillator. 47206, Cardiac pulse generator software. Class depending on the potential risk of a medical product - CLASS 3 / III. Class of electrical safety, laser safety, and other types of hazards (if any) - Not applicable for ICDs/IPGs. Safety Class of Software in accordance with IEC 62304 – C. Type of contact with the body - Long-term exposure with blood, tissues and body fluids. Sterility: sterile. Classification, depending on the possible consequences of failure during the use of the product: A.				
7. Design and appearance, photo or image, the General Scheme (flowchart)	 COMPIA MRI CRT-D DTMC2D1				
	 Torque wrench	 DF-1 Connector Pin Plug			
	8. Description of the main functional parts of a medical				
parts of a medical	Physical characteristics of implantable device: <table><tr><td>Specification</td><td>DTMC2D1</td></tr></table>			Specification	DTMC2D1
Specification	DTMC2D1				

device and technical details.

Volume	35 cm ³ ±0.5 cm ³
Mass	80 g± 1 g
H (mm) × W (mm) × D (mm)	71± 1 mm x 51± 1 mm x 13± 1 mm
Surface area of device can	57 cm ² ± 1 cm ²
Radiopaque ID*	PFZ
Radiopaque ID Medtronic*	
Materials in contact with human tissue	Titanium, polyurethane, silicone rubber, titanium dioxide
Battery	Hybrid CFx lithium/silver vanadium oxide
*Radiopaque symbol is built into each device and is visible on X-ray image. The radiopaque material is Tungsten.	

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

- the Compia MRI device
- a SureScan right atrial pacing lead or a Model 6725 pin plug for the right atrial port
- a SureScan left ventricular pacing lead
- a SureScan defibrillation lead

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Pre-implant consideration

Patient evaluation for the implant of Compia MRI different models should include the following consideration about a concomitant implant with a neurostimulator:

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procedure. For information about how to contact Medtronic, see the telephone numbers and addresses provided on the back cover of this manual.

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- CRT-D bipolar and quadripolar devices;
- SW034 Programmer Software Application;
- CareLink Programmer Model 2090;
- CareLink Encore Programmer Model 29901;
- Conexus Activator Model 27901;
- CareLink Monitor Model 2490C;
- MyCareLink Monitor Model 24950;
- 2020B CareLink Express Monitor;
- Commercially available MRI SureScan Right Atrial Pacing Leads;
- Commercially available MRI SureScan Right Ventricular Defibrillation Leads;
- Commercially available MRI SureScan Left Ventricular Pacing Leads;
- Model 6725 atrial pin plug (for device port if no right atrial pacing lead is connected);
- The same commercially available implant support instruments and accessories used with the predecessor Viva/Brava family of CRT-D device systems.

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Device name	Material Type
Digital implantable cardioverter defibrillator with cardiac resynchronization therapy and SureScan Technology COMPIA MRI CRT-D SureScan, model DTMC2D1	Left Side Shield - Titanium ASTM F67-06 Class 1 Right Side Shield - Titanium ASTM F67-06 Class 1 Fastening Clamp - Titanium ASTM F67-06 Class 1 Pin Connector Attachment Pin - Titanium Class 1 ASTM F67-06 Connector cap - Titanium class 2 F67-06 Adhesive sealant - Silicone rubber adhesive, polydimethylsiloxane MED 2000 Pin Connector Module - Polyurethane 75D, Elasthane 75D Split elastic sleeve - Silicone rubber ETR, elastomer "Silastic" class ETR Q7-4735 IS4 pin connector gasket - ZhSK (liquid silicone rubber) Silastic with pigment 2% titanium dioxide, Q7-4850
Torque wrench:	Handle - Polyetherimide G.E. Ultem 1100F-8104 Rod - Carpenter Custom 455 Stainless Steel
Pin plug Medtronic DF-1:	Handle - liquid silicone rubber (HSC) of the biomedical class "Silastic" Stem - 316 L stainless steel
Primary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated

		Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
	Secondary sterile packaging	Tray - copolymer of polyethylene terephthalate glycol (PETG) Glidex uncoated Film - Tyvek® 1073B polyethylene coated, with an adhesive layer Oliver 18B Printed Ink: -Dye Black: a. Prime Mark / b. Ultra V / c. TR4085plus™ - Dye Blue Label: Inkjet paper Adhesive: 3M P1212 Acrylic Adhesive
11. Information about Labelling and Packaging of medical device	Package content - 1 implantable cardioverter defibrillator - 1 torque wrench - 1 DF-1 pin plug - Device Manual; - Limited warranty; - Warnings and precautions related to medical procedures and EMF; - MRI Technical Manual; - Explanation of symbols for ICT and CRT-D devices; - Document "Electronic operating instructions / manuals"; - CD-ROM. Information on labelling, packaging Sterile Package: Overall dimensions, for all product variants, L x W x H, cm: 16.9 cm ± 1 cm x 13 cm ± 0.5 cm x 2.1 cm ± 0.4 cm; The mass of the package with the device, g: 120 ± 5%.	

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	Medtronic USA, Inc. 710 Medtronic Parkway N.E. Minneapolis, Minnesota 55432, USA Providing information to patients: Free call: +1 800 551 5544 (7 a.m. - 6 p.m., Monday-Friday, Central Standard Time) Fax: +1 763 514 4000		Representative office in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41 Phone: +7 495 580 7377 Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet sites: www.medtronic.eu, www.medtronic.ru Medtronic office locations worldwide are available at: http://www.medtronic.com
	Authorised representative in the European Community Medtronic B.V. Earl Bakkenstraat 10, 6422 PJ Heerlen, Netherlands +31 45 566 8000		Europe / Middle East / Africa Medtronic International Trading Sàrl Route du Molliau 31 Case Postale 84, CH-1131 Tolochenaz, Switzerland +41 21 802 7000
	Australia Medtronic Australasia Pty Ltd 97 Waterloo Road North Ryde, NSW 2113 Australia		Canada Medtronic of Canada Ltd 99 Hereford Street Brampton, Ontario L6Y 0R3 Canada +1 905 460 3800
	Waste class – B.		
15. Shelf life period	18 months.		
16. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. Store and transport the package between –18 °C and +55 °C.		
17. Requirements for maintenance, technical service and repair of a medical device	Not applicable.		
18. Warranty	Guaranteed shelf life -18 months; Warranty period - 5 years.		
19. The list of national and international regulatory documents/standards	REFERENCED STANDARDS & REVISIONS		COMPLIANCE LEVEL
	EN 556-1 Sterilization of medical devices -Requirements for medical devices to be labeled “Sterile”- Part 1: Requirements for terminally sterilized medical devices		Full Compliance
	EN 980 Symbols for use in the labelling of medical devices		Full Compliance

EN 1041 Information supplied by the manufacturer of medical devices	Full Compliance
ISO 5841-3 Implants for surgery – Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	Full Compliance
EN ISO 10993-1 Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Full Compliance
EN ISO 10993-7 Biological evaluation of medical devices - Part 7: Ethylene oxide sterilization residuals	Full Compliance
EN ISO 11135 Sterilization of health-care products - Ethylene oxide - Requirements for the development, validation and routine control of a sterilization process for medical devices	Full Compliance
ISO 11318 Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators - Dimensions and test requirement	Full Compliance
EN ISO 11607-1 Packaging for terminally sterilized medical devices – Part 1: Requirements for materials, sterile barrier systems and packaging systems	Full Compliance
EN ISO 14971 Medical devices - Application of risk management to medical devices	Full Compliance
ISO 27186 Active implantable medical devices – Four-pole connector system for implantable cardiac rhythm management devices – Dimensional and test requirements	Full Compliance
EN 45502-1 Active implantable medical devices - Part 1: General requirements for safety, marking and information to be provided by the manufacturer	Full Compliance
EN 45502-2-1 Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia (cardiac pacemakers)	Full Compliance to Section 6
EN 45502-2-2 Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Full Compliance
EN 60601-1 Medical electrical equipment -Part 1: General requirements for basic safety and essential performance	Full Compliance to Clause 14
EN 60601-1-6	Full Compliance

	Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral Standard: Usability	
	IEC 62304 Medical device software. Software life-cycle processes	Full Compliance
	EN 62366-1 Medical device-Application of usability engineering to medical devices	Full Compliance
20. Customer Service	Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic USA, Inc. Information for patients: Toll-free +1 800 551 5544 (7 a.m.– 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company in Moscow: LLC Medtronic Russia, 123112, Moscow, Presnenskaya nab., 10, floor 9, office III, room 41. Tel.: +7 495 580 7377; Fax: + 7 495 580 7378. E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com	

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