Juris Vanags 09 SEP 2019

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Date

Addendum to Device Manual

1. Name of medical product	Digital dual chamber implantable cardioverter defibrillators Evera				
	Digital dual chamber implantable cardioverter defibrillator Evera MRI XT DR SureScan, model DDMB2D1: - implantable cardioverter defibrillator; - torque wrench; - DF-1 pin plug; - Product Documentation and CD.				
2. Prescription of a medical product	The Evera MRI XT DR SureScan dual chamber implantable cardioverter defibrillator (ICD) is a multiprogrammable cardiac device that monitors and regulates the patient's heart rate by providing single or dual chamber, rate-responsive bradycardia pacing; ventricular tachyarrhythmia therapies; and atrial tachyarrhythmia therapies. The device can detect ventricular tachyarrhythmias (VT/VF) automatically and can provide treatment with defibrillation, cardioversion, and antitachycardia pacing therapies. The device can also detect atrial tachyarrhythmias (AT/AF) automatically and can provide treatment with cardioversion and antitachycardia pacing therapies. The device responds to bradyarrhythmias by providing bradycardia pacing therapy. The device also provides diagnostic and monitoring information that assists with system evaluation and patient care. The Evera dual chamber ICD devices that are the subject of this submission are listed below in: Table 1: <table border="1"> <thead> <tr> <th>Models Name</th><th>Model Number</th></tr> </thead> <tbody> <tr> <td>Evera MRI XT DR SureScan</td><td>DDMB2D1</td></tr> </tbody> </table> Application: Cardiology.	Models Name	Model Number	Evera MRI XT DR SureScan	DDMB2D1
Models Name	Model Number				
Evera MRI XT DR SureScan	DDMB2D1				

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Features/Packages	Evera MRI XT DR DDMB2D1
ICD Connectors	
1xDF4	
1xIS-1 + 1xDF4	
2xIS-1 + 2xDF-1	√
1xIS-1 + 2xDF-1	
Other Mechanical Components	
35J Defib Energy	√
Conexus® Wireless Telemetry	√
Patient Alert	√
Common System Features	
Ventricular Capture Mgmt	√
Atrial Capture Mgmt	√
Remaining Longevity Estimate	√
OptiVol® 2.0	√
AT/AF Alerts	√
AF Alerts	
Silent Alerts	√
Lead Integrity Alert	√
Brady Features, Modes and Rates	
Minimal Ventricular Pacing	√
Rate Profile Optimization	√
Rate Drop Response	√
Extended Upper Rates (UTR, USR) (to175)	√
Conducted AF Response	√

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Ventricular Arrhythmia Prevention/Suppression	
Ventricular Rate Stabilization (VRS)	√
Post-VT/VF Shock Rate and Duration	√
Ventricular Arrhythmia Detection	
VF Detection	√
VT Detection	√
FVT Detection	√
VT Monitor	√
VT Stability	√
Onset Criterion	√
PR Logic	√
High Rate Timeout	√
Wavelet	√
TWave Discrimination	√
Lead Noise Discrimination	√
VF High Rate Timeout	√
Confirmation+	√
Ventricular Arrhythmia Therapies	
Ventricular Defibrillation Therapy	√
Ventricular Cardioversion Therapy	√
Ventricular ATP Therapy	√
ATP During Charging™ + Charge Saver®	√
Atrial Intervention Pacing Features	
Post Modeswitch Overdrive Pacing	√
Atrial Rate Stabilization	√

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Medtronic, Inc.

- Confidential

Atrial Preference Pacing	√
Atrial Arrhythmia Detection	
Atrial Arrhythmia Detection (2 zone)	√
Atrial Arrhythmia Monitor (1 Zone)	√
Ventricular Based Atrial Fibrillation (VBAF)	
Atrial Arrhythmia Therapy	
Atrial Cardioversion Therapy	√
Atrial 50 Hz Burst In-Office	√
Atrial 50 Hz Burst Ambulatory	√
ATP Therapy (AT)	√
Reactive ATP [®]	√
Follow-up	
AT/AF Extended Monitoring	√
AF Extended Monitoring	
Cardiac Compass [®]	√
Hi-Fi Patient Activated EGM	√
Pre-Arrhythmia EGM Storage	√
Flashback - Atrial Episodes	√
Flashback - Ventricular Episodes	√

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Conexus wireless telemetry – This feature enables wireless transmission of data between an implanted device and the programmer in the hospital or clinic and between an implanted device and a home monitor in the patient's home.

Emergency therapies – During a patient session, defibrillation, cardioversion, fixed burst pacing, and emergency VVI can be initiated manually to treat ventricular tachyarrhythmia episodes quickly.

Live Rhythm Monitor – This window on the programmer displays ECG, Leadless ECG (LECG), Marker Channel with marker annotations, and telemetered EGM waveform traces. It also displays the patient heart rate and interval in the upper left-hand corner of the window.

Checklist – This feature presents an interactive list of common tasks that are performed during an implant session or a follow-up session. When a clinician selects a task, the associated programmer screen for that task appears. Clinicians can set up their own

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Medtronic, Inc.

- Confidential

checklists or use a Medtronic standard checklist supplied with the programmer.

Leadless ECG – This device feature allows clinicians to view and record a signal equivalent to an ECG without attaching surface ECG leads.

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OptiVol 2.0 fluid status monitoring – This feature identifies a potential increase in thoracic fluid, which may indicate lung congestion, by monitoring changes in thoracic impedance. If the change exceeds the programmed threshold, an alert tone sounds to notify the patient to seek medical attention.

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Rate Profile Optimization – This feature monitors the patient's daily and monthly sensor rate profiles and adjusts the rate response curves over time to achieve a prescribed target rate profile. The goal is to ensure that the rate response remains appropriate for the full range of patient activities.

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AT/AF detection – This feature analyzes the atrial rate and its effect on the ventricular rhythm to determine whether the patient is currently experiencing an atrial tachyarrhythmia. Evidence of an atrial tachyarrhythmia is based on the number and timing of atrial events during ventricular intervals. Depending on programming, the device delivers a programmed sequence of atrial therapies or continues monitoring without delivering therapy.

VT/VF detection – This feature uses programmable detection zones to classify ventricular events. If the number of tachyarrhythmia events in a zone exceeds a programmed threshold, the device detects a ventricular tachyarrhythmia episode. Depending on programming, the device delivers a scheduled therapy, re-evaluates the patient's heart rhythm, and terminates or re-detects the episode.

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Atrial antitachycardia pacing (ATP) – These therapies respond to an AT/AF episode or a Fast AT/AF episode with rapid sequences of pacing pulses to terminate detected atrial tachyarrhythmias. The device uses 3 programmable therapies, Burst+, Ramp,

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Programmable Active Can and SVC electrodes – The device provides the capability to disable either the Active Can or the SVC electrode as part of the high-voltage therapy delivery pathway.

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Additional operations

MRI SureScan feature – This feature allows patients to be safely scanned by an MRI machine when used according to the

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specified MRI conditions for use. Refer to the MRI technical manual for additional information.

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DDDR and DDD modes – When pacing in the DDDR or DDD mode, if the device senses an intrinsic atrial event, it schedules a ventricular paced event in response. The delay between the sensed atrial event and the corresponding ventricular paced event is the programmed Sensed AV interval. The delay between the paced atrial event and the corresponding ventricular paced event is the Paced AV interval. If the current pacing interval ends before the device senses an atrial event, the device paces the atrium and then schedules a ventricular paced event after the Paced AV interval. If a ventricular sensed event occurs during the Sensed AV interval or the Paced AV interval, ventricular pacing is inhibited. DDDR pacing occurs at the sensor rate, unless the intrinsic atrial rate is faster. DDD pacing occurs at the programmed lower rate, unless the intrinsic atrial rate is faster.

DDIR and DDI modes – When the device is pacing in the DDIR or DDI mode, sensed atrial events are not tracked. If the device senses an intrinsic atrial event, atrial pacing is inhibited, but a Sensed AV interval is not started. Instead, ventricular pacing is delivered at the current pacing rate. If the current pacing interval ends before the device senses an atrial event, the device paces the atrium and then schedules a ventricular paced event after the Paced AV interval. If a ventricular sensed event occurs during the Paced AV interval, ventricular pacing is inhibited. DDIR pacing occurs at the sensor rate. DDI pacing occurs at the programmed lower rate.

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	AAIR and AAI modes – In the AAIR and AAI modes, the atrium is paced if no intrinsic atrial events are sensed before the current pacing interval ends. Pacing occurs at the sensor rate in the AAIR mode and at the programmed lower rate in the AAI mode. AOO mode – The AOO mode provides atrial pacing at the programmed lower rate with no inhibition by intrinsic atrial events. In the AOO mode, no atrial sensing occurs. MVP mode – MVP (Managed Ventricular Pacing) is an atrial-based pacing mode that is designed to switch to a dual chamber pacing mode during AV block. MVP allows the device to function in the AAIR or AAI pacing mode when AV conduction is intact but to switch to DDDR or DDD pacing mode during AV block. It also allows the device to perform periodic conduction checks, with the ability to switch back to AAIR or AAI mode when AV conduction resumes. MVP provides backup ventricular support for transient loss of AV conduction. The AAIR<=>DDDR mode switches between the AAIR and DDDR modes and the AAI<=>DDD mode switches between the AAI and DDD modes. ODO mode – The ODO mode does not deliver ventricular or atrial pacing outputs, regardless of the intrinsic rate. The ODO mode is intended only for those situations where the clinician wants to turn off bradycardia pacing outputs from the device.	
4. Conditions for use	In medical and preventive treatment facility. The cardioverter defibrillators should be used only by or under the supervision of physicians trained.	
5. Composition and list of parts and accessories	Digital dual chamber implantable cardioverter defibrillator Evera MRI XT DR SureScan, model DDMB2D1: -1 implantable cardioverter defibrillator -1 torque wrench -1 DF-1 pin plug -Product Documentation* - CD-ROM *Product documentation includes - Device Manual, Limited warranty, Medical Procedure and EMI Precautions, Explanation of Symbols for ICD and CRT-D devices, Document: Electronic instructions for use/manuals, Appendix to the Device Manual (not in package content, printed and supplied locally on the territory of Russian Federation to the customers), Stickers, Reference manual (not in package content, available on http://manuals.medtronic.com/)	
6. Technical details of medical product	Overview of Physical Characteristics of the Evera MRI models:	
	Evera MRI Surescan XT DR	DDMB2D1
	Physical characteristics	
	Volume ³	33 cm ³ ± 0.5 cm ³
	Mass	77 g ± 1 g
	H x W x D	66 mm (±1.0 mm) x 51 mm (±1.0 mm) x 13 mm (±0.5 mm)

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Surface area of device can	57 cm ² ± 0.5 cm ²
Radiopaque ID ^b	PFZ*
Medtronic radiopaque identifier ^b	
Materials in contact with human tissue ^c	Titanium, polyurethane, silicone rubber
Battery	Hybrid CFx lithium/silver vanadium oxide

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


- | | |
|--|--|
| | <ul style="list-style-type: none">1 DF-1 connector port, SVC2 DF-1 connector port, RV3 Device Active Can electrode4 IS-1 connector port, RV5 IS-1 connector port, A6 Suture holes |
|--|--|

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Battery characteristics

Manufacturer	Medtronic Energy and Component Center
Model	M945899A
Chemistry	Hybrid CFx lithium/silver vanadium oxide

Battery electrical specifications

Nominal voltage	3.2 V
Mean capacity to RRT	1.0 Ah
Minimum capacity after RRT	0.1 Ah

Device electrical specifications

Pacing rate limit (protective feature) ^a	200 min ⁻¹
Input impedance	150 kΩ minimum

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
DF-1 Connector Pin Plug (unit: mm)

DF-1 Connector Pin Plug included in package with the following models only: Evera MRI DDMB2D1. This device is a plug to fill an IS-1 (ISO 5841-3) compliant connector port of pulse generators, such as an implantable cardiac pacemaker, to which an implantable IPG/ICD lead connector is not connected.

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Medtronic, Inc.

- Confidential

	Model 27901 Conexus Activator – The Medtronic Model 27901 Conexus Activator allows you to turn on Conexus wireless telemetry for implanted devices that support wireless telemetry. The Conexus Activator is used with the Medtronic Analyzer Programmer with Conexus telemetry in the hospital or clinic. Medtronic Analyzer – The system supports the use of the Medtronic Analyzer, an accessory of the Medtronic CareLink Programmer. The system allows you to have a device session and an analyzer session running at the same time, to switch quickly from one session to the other without having to end or restart sessions, and to send data from the analyzer to the programmer. Medtronic Patient Monitor – Patients use the Medtronic patient monitor to gather information automatically from their implanted devices and communicate the information to their physicians through the Medtronic CareLink Network. For information on using the patient monitor, refer to the patient manual; for connection and usage information, refer to the monitor literature. CareLink Network – The Medtronic CareLink Network is an Internet-based service that allows a patient to transmit full cardiac device information from home or other locations to the physician over a secure server. The CareLink Network may be unavailable in some geographic locations. CareLink Express monitor Model 2020B – The Medtronic CareLink Express monitor Model 2020B is a portable interrogation and data transfer tool used with Medtronic implantable cardiac devices. Model 2696 InCheck Patient Assistant – Patients can use the Medtronic Model 2696 InCheck Patient Assistant to initiate recording of cardiac event data in the device memory, verify whether the implanted device has detected a suspected atrial tachyarrhythmia, and, if the device is programmed to allow patient-activated cardioversion, request delivery of atrial cardioversion therapy. Model 2490C CareLink Monitor – Patients use the Medtronic Model 2490C CareLink Monitor to gather information automatically from their implanted devices and communicate the information to their physicians. For information on using the CareLink Monitor, refer to the patient manual; for connection and usage information, refer to the monitor literature.																		
8. EMC information	<table><tr><th colspan="3">Electromagnetic emission</th></tr><tr><td colspan="3">Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment</td></tr><tr><td>Electromagnetic emission test</td><td>Conformity</td><td>Electromagnetic environment – instructions</td></tr><tr><td>Radio interference according to CISPR 11</td><td>Group 1</td><td rowspan="4">Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.</td></tr><tr><td>Radio interference to CISPR 11</td><td>Class B</td></tr><tr><td>Harmonic components of the consumed current in accordance with IEC 61000-3-2</td><td>Do not use</td></tr><tr><td>Voltage fluctuations and flicker according to IEC 61000-3-3</td><td>Do not use</td></tr></table>	Electromagnetic emission			Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment			Electromagnetic emission test	Conformity	Electromagnetic environment – instructions	Radio interference according to CISPR 11	Group 1	Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.	Radio interference to CISPR 11	Class B	Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use	Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use
Electromagnetic emission																			
Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment																			
Electromagnetic emission test	Conformity	Electromagnetic environment – instructions																	
Radio interference according to CISPR 11	Group 1	Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.																	
Radio interference to CISPR 11	Class B																		
Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use																		
Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use																		

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Noise immunity			
Immunity test	The test level of IEC 60601	Level of compliance	Electromagnetic environment - instructions
Electrostatic Discharge (ESD) in accordance with IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%
Magnetic field of industrial frequency according to IEC 61000-4-8	3 A/m	0,3 A/m	The magnetic fields of the industrial frequency must be measured at the designated installation site to ensure that the field strength is sufficiently low
Immunity test	The test level of IEC 60601	The actual level of noise immunity	Level of compliance
Conductive noise induced by radio-frequency electromagnetic fields in accordance with IEC 61000-4-6	3 V (rms) in the band from 150 kHz to 80 MHz outside the frequencies allocated for the LPM of HF devices	1.7 V (rms value)	1 V (rms value)
The emitted radio-frequency electromagnetic field according to IEC 61000-4-3	10 V (RMS) in the band from 150 kHz to 80 MHz at the frequencies allocated for the LPM of high-frequency devices	2.3 V (rms value)	1 V (rms value)
		13 V/m	10 V/m
Noise immunity			
The test for Noise immunity	Test Level according to IEC 60601	Level Conformity	Electromagnetic environment - instructions
Electrostatic Discharge (ESR) according to IEC 61000-4-2	± 6 kV – contact discharge ± 8 kV – air discharge	± 6 kV – contact discharge ± 8 kV – air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%.

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	A magnetic field Industrial Frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of the industrial frequency must correspond to the typical conditions of a commercial or hospital environment.
	Radio frequency Electromagnetic field IEC 61000-4-3	3 V/m In the band from 80 MHz to 2.5 GHz	3 V/m	However, the readings of the weights built into the bed do not have a critical value for the operation of a medical device: «Digital Implantable dual chamber Cardioverter Defibrillators Evera» $D = 1,2\sqrt{P}$ in the band from 26 to 800 MHz $D = 2,3\sqrt{P}$ in the band from 800 MHz to 2.5 GHz, where P is the nominal maximum output transmitter power in Watts (W) set by the manufacturer, and d - Recommended spatial clearance in meters (m). The field strength in the propagation of radio waves from stationary radio transmitters from the results of observations of the electromagnetic situation should be lower than the level of compliance in each frequency band. Interference may occur near equipment  marked with a sign. This device may not cause harmful interference, and this. To accept any interference received, including interference that may cause undesired operation of the device. This device must not cause harmful interference to radio stations operating in the 400.150-406.000 MHz band and used for meteorological purposes, for meteorological satellites and Systems of Earth exploration satellites. This device must accept any interference received, including interference that may cause undesired operation.

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Material	Device Components
Titanium grade 1 <u>ASTM F67</u>	LH shield
Fe to 0.12	RH shield
C to 0.05	Connect attach pin
Si to 0.08	Fastener bracket
N to 0.03	
Ti 99.5 – 99.99	
O to 0.1	
H to 0.008	
Titanium grade 2 <u>ASTM F67</u>	End cap
Fe to 0.18	
C to 0.07	
Si to 0.1	
N to 0.03	
Ti 98.61 – 99.7	
O to 0.12	
H to 0.01	
Silicone rubber MED 2000	Adhesive
Polyurethane 75D Elasthane 75D [™] /Tecothane TT-1075D-M	Connector module
Silicone rubber ETR Elastomer Q7-4735	Split grommet
LSR with titanium dioxide 2% (Q74850, colorant Ferro white 61-182000)	IS4 seal

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
Torque wrench and DF-1 pin plug (for Evera XT DR MRI DDMB2D1):

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11. Information on labelling, packaging

Products supplied in packaging marked with the following information:

- product name and catalog number;
- LOT of the product;
- trademark and address of the manufacturer;
- single-use;
- sterile, do not resterilize;
- label complies with European Council 90/385/EEC;
- manufacturing date and shelf life of the product;
- storage conditions;
- do not use if the package is damaged;
- characteristics of the product

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2000-00-00
YYYY-MM-DD

www.medtronic.com/manuals

Medtronic, Inc., 710 Medtronic Parkway
Minneapolis, MN 55432 USA
Manufactured in
Tessenderlo, Switzerland
www.medtronic.com

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
Medtronic

DF-1 Evera MRI™ XT DR SureScan™ DDMB2D1 DDE-DDDR

Adaptive	AAI ↔ DDD	1 3.8 V @ → 0.4 ms @	36/42 J
	↑ 130 min ⁻¹ ↓ 60 min ⁻¹ Upper tracing rate and lower rate	Appropriate and pulse with RA	Therapy delivered in 30 seconds VF
	0.3 mV Detecting RA	1 3.8 V @ → 0.4 ms @	ON Active Can™
	0.3 mV Detecting RV	Appropriate and pulse with RV	OFF
	→ 150 ms → 150 ms A/V interval (paced and sensed)	OFF Detecting VT, PVT, VF	Therapy AT A ²
		OFF Detecting (RV) VT, PVT	OFF Therapy AT A ²

35°C / 95°F to 131°F / 268°F
 FCC Class II LFPD COMPLIANT (IC 2400-AB-01-0013)
STERILE **EO**
 Medtronic, Inc. 710 Medtronic Parkway
 Minneapolis, MN 55432, USA
 CE
 Medtronic B.V. Eerste Bakkensloot 10
 1142 ZJ Haren, The Netherlands
www.medtronic.com
 Manufactured in
 Tuzlukent, Switzerland
 Product of Switzerland
CE0123
 Medtronic A - Medtronic A

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


	Use by
	Do not use if package is damaged
	Date of manufacture
	Do not reuse
	Authorized representative in the European Community
AAI ⇔ DDD	AAI ⇔ DDD mode
	Package contents
	Torque wrench
	Pin plug, DF-1
	Product documentation

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
	Therapies (delivered and stored): VF
	Active Can
	Defibrillation
	AV cardioversion
	Ramp: RA, RV
	Ramp+: RV
	Burst: RV
	Burst+: RA
	50 Hz Burst
	SureScan symbol

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- Confidential

Thickness: 178 microns

4. Compatibility with the methods of sterilization with which the use of the material is intended:
Methods of sterilization:

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5. Compatibility with the processes of forming and gluing:
Termoskleivaetsya with a large number of thermoplastics, as PVC, styrene, PETG, HIPS; With polyethylene some kinds of polypropylene.

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12. pH: Tyvek has a neutral pH of 7.

13. Adhesive characteristics:
Peel strength (strength of the adhesive layer): 4.53 H-18.14 H
Weight of the adhesive layer: $4.5 \text{ kg} \pm 0.7 \text{ kg} / 279 \text{ m}^2$
Adhesion transfer: > 75%
Shelf life: <48 months
White color
The density of glue application: > $1.03 \text{ (g / cm}^3\text{)}$

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
18. Used dyes:
Inks are black C70, produced by Sun Chemical Corp., USA;
Blue SUNSPERSE 15, produced by Sun Chemical Corp, USA

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
Packaging Description
Shelf Box – Secondary Packaging

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
Package Part	Material
Sterile tray	Pre-shaped (thermo-formed) uncoated PETG (polyethylene terephthalate, glycol modified)
Lid	Spun bonded olefin lids coated with Oliver 18B adhesive

* The deviation of the technical characteristics of the package is $\pm 5\%$, unless otherwise specified.

12. Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)

Sterilization&Validation	
Sterilization Method	Ethylene Oxide (EO)
Sterilization cycle/Process name	30-minute 100% EO sterilization process
Applicable Sterilization Process Validation Standard	ISO 11135-1:2007
Sterilant Gas type	100% EO
Relative Humidity (RH)	55.0-85.0%
EO Gas Concentration	≥ 720 mg/l
Temperature during EO Gas exposure	45.0-55.0 °C
Pressure during EO gas Exposure	55.0-75.0 kPa
EO Gas Exposure time	27-33 min

Medtronic has sterilized the package contents with ethylene oxide before shipment (30 min. 100 % EO). A SAL in excess of 10^{-6} is achieved. Detected EO residual levels were less than 4 mg and detected ECH residual levels were less than 9 mg. Therefore, compliance with ISO 10993-7.

13. Requirements for the protection of the environment in the application of medical products

Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See <http://recycling.medtronic.com> for instructions on proper disposal of this product.

Explant and disposal

Consider the following information related to device explant and disposal:

- To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device.
- Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures.
- Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices.
- Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. See the back cover for addresses.

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disinfection and sterilization cleaning																																			
17. Shelf life period	18 months Projected service life The projected service life in years for the device is shown in the tables below. The data is based on pacing outputs programmed to the specified mode, the specified pacing amplitude, 0.4 ms pulse width, and 60 min⁻¹ pacing rate. The service life projections are based on the following assumptions: • Semi-annual maximum energy charging frequency • Pre-arrhythmia EGM storage programmed to On for a 6-month period (two 3-month follow-up intervals), over the entire life of the device • 3 hours of wireless telemetry during implant • A quarterly schedule of Medtronic CareLink Monitor remote transmissions • 1 hour of in-office wireless telemetry annually • Typical shelf storage time before implant Projected service life estimates are based on accelerated battery discharge data and device modeling as specified. Do not interpret these values as precise numbers. Projected service life in years Projected service life in years (Evera MRI XT DR SureScan DDMB2D1) <table><tr><th rowspan="2">Pacing mode, percent pacing</th><th rowspan="2">Pacing amplitude</th><th colspan="2">Projected service life in years</th></tr><tr><th>Pacing impedance 500 Ω</th><th>Pacing impedance 600 Ω</th></tr><tr><td rowspan="2">DDD, 0%</td><td>2.5 V</td><td>9.7</td><td>9.7</td></tr><tr><td>3.5 V</td><td>9.6</td><td>9.6</td></tr><tr><td rowspan="2">DDD, 15%</td><td>2.5 V</td><td>9.1</td><td>9.2</td></tr><tr><td>3.5 V</td><td>8.7</td><td>8.8</td></tr><tr><td rowspan="2">DDD, 50%</td><td>2.5 V</td><td>8.3</td><td>8.5</td></tr><tr><td>3.5 V</td><td>7.1</td><td>7.4</td></tr><tr><td rowspan="2">DDD, 100%</td><td>2.5 V</td><td>7.3</td><td>7.6</td></tr><tr><td>3.5 V</td><td>5.7</td><td>6.1</td></tr></table>	Pacing mode, percent pacing	Pacing amplitude	Projected service life in years		Pacing impedance 500 Ω	Pacing impedance 600 Ω	DDD, 0%	2.5 V	9.7	9.7	3.5 V	9.6	9.6	DDD, 15%	2.5 V	9.1	9.2	3.5 V	8.7	8.8	DDD, 50%	2.5 V	8.3	8.5	3.5 V	7.1	7.4	DDD, 100%	2.5 V	7.3	7.6	3.5 V	5.7	6.1
Pacing mode, percent pacing	Pacing amplitude			Projected service life in years																															
		Pacing impedance 500 Ω	Pacing impedance 600 Ω																																
DDD, 0%	2.5 V	9.7	9.7																																
	3.5 V	9.6	9.6																																
DDD, 15%	2.5 V	9.1	9.2																																
	3.5 V	8.7	8.8																																
DDD, 50%	2.5 V	8.3	8.5																																
	3.5 V	7.1	7.4																																
DDD, 100%	2.5 V	7.3	7.6																																
	3.5 V	5.7	6.1																																
18. Operating conditions	Used in medical and preventive treatment facility. Should be used only by or under the supervision of physicians trained. internal body environment 32 °C-42 °C, humidity 80-100%.																																		
19. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. There are no other requirements for storage and transportation, except for temperature: Store and transport the package between -18 °C and +55 °C.																																		

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a) the cost of the Device has been charged to the account of the relevant Healthcare Institution; and

b) the Healthcare Institution's claim under this Limited Warranty is solely for the purpose of obtaining a second Medtronic device to be used by it as a Replacement (as defined below) for the Device for the same patient (this category of Purchaser is referred to in this Limited Warranty as an "Institution-Purchaser").

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a) Issue a credit to the Purchaser for a Medtronic device to be used as a replacement of the Device (the "Replacement") for the same patient, equal to the lesser of the net price paid by the Purchaser: (i) for the Device (the "Original Purchase Price"), or (ii) for the Replacement (the "Actual Purchase Price"); or,

b) Without charge, provide to the Purchaser a Replacement for use in the patient in whom the Device was originally implanted that is functionally comparable to the Device.

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(1) The Device must be implanted on or before its "Use By" or "Use Before" date in conjunction with Medtronic leads or leads of equal quality and comparable electrical characteristics.

(2) The Device must be used in accordance with the labeling and not altered or subjected to misuse, reuse, abuse, accident, or improper handling.

(3) In jurisdictions where it is a legal requirement or established practice to register a Device with the manufacturer, all Device

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ISO 14708-6:2010	Implants for surgery - Active implantable medical devices - Part 6: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (including implantable defibrillators)	—	Full Compliance
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices (ISO 14971:2007, corrected version 2007-10-01)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
ISO 11607-1:2009	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2006)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
ISO 11607-2:2006	Packaging for terminally sterilized medical devices Part 2: Validation requirements for forming, sealing and assembly processes - First edition	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN 1041:2008+2013	Information supplied by the manufacturer with medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 11135-1:2014	Sterilization of medical devices – Validation and routine control of Ethylene Oxide sterilization	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN 980:2008	Graphical symbols for use in the labelling of medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-1:2009/AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Harmonized under the AIMD and MD Directives of Europe	Full compliance for applicable sections: (Meets all requirements for implanted devices in tissue/bone listed in tables 1 and 2)
EN ISO 10993-2:2006	Biological evaluation of medical devices - Part 2: Animal welfare requirements		Full compliance

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EN ISO 10993-3:2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-4:2017	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood	Harmonized under the AIMD and MD Directives of Europe	Not available
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-6:2016	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-7:2008/AC:2009	Biological evaluation of medical devices - Part 7: Ethylene Oxide sterilization residuals	Harmonized under the AIMD and MD Directives of Europe	Full Compliance for applicable sections (section 4.3.2 Permanent contact device)
EN ISO 10993-9:2009	Biological evaluation of medical devices - Part 9: Framework for identification and quantification of potential degradation products	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization		Full Compliance
EN ISO 10993-11:2006	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-13:2010	Biological evaluation of medical devices - Part 13: Identification and quantification of degradation products from polymeric medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-15:2000	Biological evaluation of medical devices. Part 15: Identification and quantification of degradation products from metals and alloys	Harmonized under the AIMD and MD Directives of Europe	Full Compliance

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EN ISO 10993-16:2010	Biological evaluation of medical devices - Part 16: Toxicokinetic study design for degradation products and leachables	Harmonized under the AIMD and MD Directives of Europe	The evaluation was not carried out due to the availability of data on the use of these materials in marketable products and biosafety testing
EN ISO 10993-17:2002	Biological evaluation of medical devices, Part 17. Establishment of allowable limits for leachable substances	Harmonized under the AIMD and MD Directives of Europe	The evaluation was not carried out due to the availability of data on the use of these materials in marketable products and biosafety testing
EN 556-1: 2001/AC:2006	Sterilization of medical devices - Requirements for medical devices to be labeled "Sterile" - Part 1: requirements for terminally sterilized medical devices	—	Full Compliance
ISO 5841-3: 2013	Implants for surgery - Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	—	Full compliance for applicable sections (Section 4.2 is applicable to leads only and does not apply)
ISO 11318	Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators – Dimensions and test requirements	—	Full Compliance
EN 62366-1:2015	Medical devices - Application of usability engineering to medical devices	Harmonized under the MD Directives of Europe	Full compliance only for user interface and other design outputs developed on the Evera MRI project
EN 62304:2016	Medical device software. Software life-cycle processes	Harmonized under the AIMD and MD Directives of Europe	Full compliance only for code developed on the Evera MRI project
EN 60601-1: 2006 AC:2010	General requirements for safety of Medical Electrical devices	Harmonized under the AIMD and MD Directives of Europe (1st Edition)	Clause 14 only (programmer software)

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↑ 3.5 V @
 → 0.4 ms @
 Asynchronous and pulse width RA
 ↑ 3.5 V @
 → 0.4 ms @
 Asynchronous and pulse width RV
 OFF
 Detect for VT, VVI, W
 OFF
 Detect for RVO, VVI, FVT

3542 J
 Therapy (synchronous and asynchronous) W
 ON
 Active Can™
 OFF
 Therapy AT, AF
 OFF
 Therapy AT, AF

55°C
131°F
 18°C
0°F

FDC IF: LPS/MS/IMPLANTS
 IC: 3408D IMPLANTS

STERILE EO

SN WWW000000W
 2000-00-00
 YYY-AAA-DD
 2000-00-00
 YYY-AAA-DD

www.medtronic.com/manuals


Medtronic, Inc., 710 Medtronic Parkway
 Minneapolis, MN 55432 USA
 Manufactured in:
 Touchez, Switzerland
 www.medtronic.com

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The diagram illustrates the components of the Active Can™ system and their interconnections. On the left is a large oval representing the 'Can' with a lightning bolt symbol inside a triangle. To its right is a circular 'ACTIVE CAN™' unit with three internal circles labeled 'A', 'BVC', and 'RV'. Above this unit are three small circular icons: a 'Sunlight' icon, a 'Rain' icon, and an 'ICD (Flow + Handberry)' icon. Below the 'ACTIVE CAN™' unit is a 'Pie Logger' icon. To the right of the 'Pie Logger' are four icons: a '1x' icon, a 'Pie Logger' icon, a 'Pie Logger' icon, and a 'Pie Logger' icon. A legend at the bottom left shows two circles: a solid circle labeled '10-1 BI' and an open circle labeled 'GP-1'. Arrows indicate the flow of information or data between these components.

De Bouveret AV-septum Hirsch, RV Hirsch, RA Hirsch, RA, RV Hirsch, RV Hirsch, RA, RV

A-V interval (paced and sensed)

Theoretical (RV) VI, FV

Intervention A

© 2004 Blackwell Publishing Ltd *Journal of Internal Medicine* 255: 111–117

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
Cindi A. Bigelbach
9/9/19

Juris Vanags 09 SEP 2019

Juris Vanags
Sr. Regulatory Affairs Specialist
Medtronic, Inc.

Date

Addendum to Device Manual

1. Name of medical product	Digital dual chamber implantable cardioverter defibrillators Evera				
	Digital dual chamber implantable cardioverter defibrillator Evera MRI S DR SureScan, model DDMC3D1: - implantable cardioverter defibrillator; - torque wrench; - DF-1 pin plug; - Product Documentation and CD.				
2. Prescription of a medical product	The Evera MRI S DR SureScan dual chamber implantable cardioverter defibrillator (ICD) is a multiprogrammable cardiac device that monitors and regulates the patient's heart rate by providing single or dual chamber, rate-responsive bradycardia pacing; ventricular tachyarrhythmia therapies; and atrial tachyarrhythmia therapies. The device can detect ventricular tachyarrhythmias (VT/VF) automatically and can provide treatment with defibrillation, cardioversion, and antitachycardia pacing therapies. The device can also detect atrial tachyarrhythmias (AT/AF) automatically and can provide treatment with cardioversion and antitachycardia pacing therapies. The device responds to bradyarrhythmias by providing bradycardia pacing therapy. The device also provides diagnostic and monitoring information that assists with system evaluation and patient care. The Evera dual chamber ICD devices that are the subject of this submission are listed below in: Table 1: <table border="1"> <thead> <tr> <th>Models Name</th><th>Model Number</th></tr> </thead> <tbody> <tr> <td>Evera MRI S DR SureScan</td><td>DDMC3D1</td></tr> </tbody> </table> Application: Cardiology.	Models Name	Model Number	Evera MRI S DR SureScan	DDMC3D1
Models Name	Model Number				
Evera MRI S DR SureScan	DDMC3D1				

3. Description of the principles underlying the work of the medical products and their features

Feature Comparison of Evera MRI XT and Evera MRI S ICD devices:
[Key: √ = feature present in device; empty box = feature absent from device]

Features/Packages	Evera MRI S DR DDMC3D1
ICD Connectors	
1xDF4	
1xIS-1 + 1xDF4	
2xIS-1 + 2xDF-1	√
1xIS-1 + 2xDF-1	
Other Mechanical Components	
35J Defib Energy	√
Conexus [®] Wireless Telemetry	√
Patient Alert	√
Common System Features	
Ventricular Capture Mgmt	√
Atrial Capture Mgmt	√
Remaining Longevity Estimate	√
AT/AF Alerts	√
AF Alerts	
Silent Alerts	√
Lead Integrity Alert	√
Brady Features, Modes and Rates	
Minimal Ventricular Pacing	√
Rate Profile Optimization	√
Rate Drop Response	√
Extended Upper Rates (UTR, USR) (to 175)	√
Conducted AF Response	√
Ventricular Arrhythmia Prevention/Suppression	

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Ventricular Rate Stabilization (VRS)	√
Post-VT/VF Shock Rate and Duration	√
Ventricular Arrhythmia Detection	
VF Detection	√
VT Detection	√
FVT Detection	√
VT Monitor	√
VT Stability	√
Onset Criterion	√
PR Logic	√
High Rate Timeout	√
Wavelet	√
TWave Discrimination	√
Lead Noise Discrimination	√
VF High Rate Timeout	√
Confirmation+	√
Ventricular Arrhythmia Therapies	
Ventricular Defibrillation Therapy	√
Ventricular Cardioversion Therapy	√
Ventricular ATP Therapy	√
ATP During Charging™ + Charge Saver®	√
Atrial Intervention Pacing Features	
Post Modeswitch Overdrive Pacing	√
Atrial Rate Stabilization	√
Atrial Preference Pacing	√
Atrial Arrhythmia Detection	

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Atrial Arrhythmia Detection (2 zone)	√
Atrial Arrhythmia Monitor (1 Zone)	√
Ventricular Based Atrial Fibrillation (VBAF)	
Atrial Arrhythmia Therapy	
Atrial Cardioversion Therapy	√
Atrial 50 Hz Burst In-Office	√
Atrial 50 Hz Burst Ambulatory	√
ATP Therapy (AT)	√
Reactive ATP [®]	√
Follow-up	
AT/AF Extended Monitoring	√
AF Extended Monitoring	
Cardiac Compass [®]	√
Hi-Fi Patient Activated EGM	√
Pre-Arrhythmia EGM Storage	√
Flashback - Atrial Episodes	√
Flashback - Ventricular Episodes	√
Programmer software features For more information about these features, see the reference manual. Conexus wireless telemetry – This feature enables wireless transmission of data between an implanted device and the programmer in the hospital or clinic and between an implanted device and a home monitor in the patient's home. Emergency therapies – During a patient session, defibrillation, cardioversion, fixed burst pacing, and emergency VVI can be initiated manually to treat ventricular tachyarrhythmia episodes quickly. Live Rhythm Monitor – This window on the programmer displays ECG, Leadless ECG (LECG), Marker Channel with marker annotations, and telemetered EGM waveform traces. It also displays the patient heart rate and interval in the upper left-hand corner of the window. Checklist – This feature presents an interactive list of common tasks that are performed during an implant session or a follow-up session. When a clinician selects a task, the associated programmer screen for that task appears. Clinicians can set up their own checklists or use a Medtronic standard checklist supplied with the programmer. Leadless ECG – This device feature allows clinicians to view and record a signal equivalent to an ECG without attaching surface ECG leads.	

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

TherapyGuide – This feature suggests a set of parameters based on the programmed information about the patient's clinical conditions. TherapyGuide does not replace a physician's expert judgment. The physician can accept, reject, or modify any of the suggested parameter values.

Patient Information – This feature allows clinicians to store patient-related information on the programmer that they can view and print during a patient session.

Diagnostic data features

When MRI SureScan is programmed to On, diagnostic data is not collected. Before performing an MRI scan, refer to the MRI technical manual for MRI-specific warnings and precautions.

For more information about these features, see the reference manual.

Quick Look II – This screen on the programmer presents overview data about device operation and patient rhythms collected since the last patient session. It includes links to more detailed status and diagnostic information stored in the device, such as arrhythmia episodes and therapies provided.

Medtronic CareAlert – If the device identifies any CareAlert programmed or automatic alert conditions, this feature sends a wireless alert signal to the CareLink Monitor, transmits an alert notification to the clinic, and sounds a patient alert tone to notify the patient to seek medical attention.

RV Lead Integrity Alert – This feature sounds an alert tone to warn the patient that a potential RV lead problem is suspected, which could indicate a lead fracture. When the alert criteria are met, device settings are automatically adjusted to avoid the delivery of inappropriate therapy.

Cardiac Compass Trends – This feature presents an overview of the patient's condition over the past 14 months with graphs that display long-term clinical trends in heart rhythm and device status, such as frequency of arrhythmias, heart rates, and device therapies.

Heart Failure Management Report – This printed report summarizes the patient's clinical status and observations since the last follow-up appointment. The report provides graphs that show trends in heart rates, arrhythmias, and fluid accumulation indicators over the past 14 months.

Arrhythmia episode data – The system compiles an arrhythmia episode log that the clinician can use to view summary and detailed diagnostic data quickly, including stored EGM, for the selected arrhythmia episode. Also available on the programmer are episode and therapy counters, stored data showing the number of times that arrhythmias and therapies have occurred.

Flashback Memory – This diagnostic feature records the intervals that immediately preceded tachyarrhythmia episodes or that preceded the last interrogation of the device and plots the interval data over time.

Rate Drop Response episodes data – This feature displays beat-to-beat data that is useful in analyzing Rate Drop Response episodes and the events leading up to these episodes. The feature records data about episodes that meet the programmed rate drop detection criteria.

Rate Histograms – This diagnostic feature shows range distributions for the patient's heart rate.

Pacing features

For more information about these features, see the reference manual.

Auto-adjusting sensitivity – This feature automatically adjusts the sensitivity thresholds after specific paced events and sensed events occur.

Rate Response – This feature adjusts the cardiac pacing rate in response to changes in sensed patient activity.

Rate Profile Optimization – This feature monitors the patient's daily and monthly sensor rate profiles and adjusts the rate response

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

curves over time to achieve a prescribed target rate profile. The goal is to ensure that the rate response remains appropriate for the full range of patient activities.

Capture Management – This feature monitors pacing thresholds with daily pacing threshold searches and, if programmed to do so, adjusts the pacing amplitudes toward a target amplitude.

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Auto PVARP – This feature adjusts the Post-Ventricular Atrial Refractory Period (PVARP) in response to changes in the patient's heart rate or pacing rate. PVARP is longer at lower tracking rates to prevent pacemaker-mediated tachycardia (PMT) and shorter at higher rates to maintain 1:1 tracking.

Rate Drop Response – This feature monitors the heart for a significant drop in rate and responds by pacing the heart at an elevated rate for a programmed duration.

Sleep – This feature causes the device to pace at a slower rate during a programmed sleep period.

Non-Competitive Atrial Pacing (NCAP) – This feature prohibits atrial pacing within a programmable interval after a refractory atrial event.

Pacemaker-mediated tachycardia (PMT) Intervention – This feature automatically detects and interrupts device-defined PMTs.

PVC Response – This feature extends PVARP following a premature ventricular contraction (PVC) to avoid tracking a retrograde P-wave and to prevent retrograde conduction from inhibiting an atrial pace.

Ventricular Safety Pacing (VSP) – This feature prevents inappropriate inhibition of ventricular pacing caused by crosstalk or ventricular oversensing.

Mode Switch – This feature switches the device pacing mode from a dual-chamber atrial tracking mode to a nontracking mode during an atrial tachyarrhythmia. This feature prevents rapid ventricular pacing that can result from tracking a high atrial rate and restores the programmed pacing mode when the atrial tachyarrhythmia ends.

Conducted AF Response – This feature helps to promote a regular ventricular rate during conducted AT/AF episodes by increasing the pacing rate in concert with the patient's intrinsic ventricular response.

Atrial Rate Stabilization (ARS) – This feature adapts the atrial pacing rate in response to a PAC (premature atrial contraction) to avoid long sinus pauses following short atrial intervals.

Atrial Preference Pacing (APP) – This feature maintains a consistent activation sequence by providing continuous pacing that is closely matched to the intrinsic rate.

Post Mode Switch Overdrive Pacing (PMOP) – This feature works with the Mode Switch feature to deliver overdrive atrial pacing during the vulnerable phase following an AT/AF episode termination.

Post VT/VF Shock Pacing – This feature provides temporary overdrive pacing for a programmed duration after a ventricular high-voltage therapy.

Ventricular Rate Stabilization (VRS) – This ventricular rhythm management feature adjusts the pacing rate dynamically to eliminate the long pause that typically follows a premature ventricular contraction (PVC).

Tachyarrhythmia detection features

When MRI SureScan is programmed to On, tachyarrhythmia detection and therapies are suspended. Before performing an MRI scan, refer to the MRI technical manual for MRI-specific warnings and precautions.

For more information about these features, see the reference manual.

AT/AF detection – This feature analyzes the atrial rate and its effect on the ventricular rhythm to determine whether the patient is

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

currently experiencing an atrial tachyarrhythmia. Evidence of an atrial tachyarrhythmia is based on the number and timing of atrial events during ventricular intervals. Depending on programming, the device delivers a programmed sequence of atrial therapies or continues monitoring without delivering therapy.

VT/VF detection – This feature uses programmable detection zones to classify ventricular events. If the number of tachyarrhythmia events in a zone exceeds a programmed threshold, the device detects a ventricular tachyarrhythmia episode. Depending on programming, the device delivers a scheduled therapy, re-evaluates the patient's heart rhythm, and terminates or redetects the episode.

PR Logic – This set of features uses pattern and rate analysis to discriminate between supraventricular tachycardias (SVTs) and true ventricular tachyarrhythmias and to withhold inappropriate VT/VF detection and therapy during episodes of rapidly conducted SVT.

Wavelet – This feature is designed to prevent detection of rapidly conducted SVTs as ventricular tachyarrhythmias by comparing the shape of each QRS complex during a fast ventricular rate to a template. The feature offers the option to collect and maintain the stored template automatically.

Onset – This feature helps prevent detection of sinus tachycardia as VT by evaluating the acceleration of the ventricular rate.

Stability – This feature helps to prevent detection of atrial fibrillation as ventricular tachyarrhythmia by evaluating the stability of the ventricular rate. If the device determines that the ventricular rate is not stable, it withholds VT detection.

High Rate Timeout – This feature allows the device to deliver therapy for any ventricular tachyarrhythmia that continues beyond the programmed length of time.

TWave Discrimination – This feature withholds VT/VF detection when a fast ventricular rate is detected because of oversensed T-waves, avoiding the delivery of an inappropriate therapy.

RV Lead Noise Discrimination – When the device identifies lead noise due to a suspected lead problem, this feature withholds VT/VF detection to prevent the delivery of inappropriate therapy. Also, if programmed, an alert tone sounds to notify the patient to seek medical attention.

Tachyarrhythmia therapy features

When MRI SureScan is programmed to On, tachyarrhythmia detection and therapies are suspended. Before performing an MRI scan, refer to the MRI technical manual for MRI-specific warnings and precautions.

For more information about these features, see the reference manual.

Atrial therapy scheduling – This feature enables the clinician to program the delivery of automatic atrial therapies. Each time that an AT/AF therapy is needed, the device schedules one of the available therapies based on clinician programming. The Reactive ATP option allows the device to repeat programmed atrial antitachycardia pacing (ATP) therapies during long AT/AF episodes. Therapies are repeated after a programmed interval or when the atrial rhythm changes in regularity or cycle length.

Atrial antitachycardia pacing (ATP) – These therapies respond to an AT/AF episode or a Fast AT/AF episode with rapid sequences of pacing pulses to terminate detected atrial tachyarrhythmias. The device uses 3 programmable therapies, Burst+, Ramp, and 50 Hz Burst, to treat an episode.

Atrial cardioversion – This therapy delivers a high-voltage shock to treat an AT/AF episode or a Fast AT/AF episode. Atrial cardioversion delivery is synchronized to a sensed ventricular event and cannot exceed a programmable daily limit within programmable times. The patient also can request atrial cardioversion using an external Patient Assistant. Patient-activated atrial cardioversion is delivered only if an AT/AF episode is detected at the time of the request.

Programmable Active Can and SVC electrodes – The device provides the capability to disable either the Active Can or the SVC

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

electrode as part of the high-voltage therapy delivery pathway.

Ventricular fibrillation (VF) therapies – Automatic defibrillation shocks are available to treat VF episodes. The first defibrillation therapy requires VF confirmation before delivery. After the first shock has been delivered, shocks are delivered asynchronously if synchronization fails. ATP During Charging allows the device to deliver a ventricular antitachycardia pacing therapy sequence while the device charges its capacitors for the first defibrillation therapy. The device can also be programmed to attempt an additional ATP therapy sequence before charging starts.

Ventricular antitachycardia pacing (ATP) – These therapies respond to a VT episode or an FVT episode with rapid sequences of pacing pulses to terminate detected ventricular tachyarrhythmias. Therapy options include Burst, Ramp, and Ramp+, each with a programmable number of sequences.

Ventricular cardioversion – This therapy delivers a high-voltage shock to treat a VT or FVT episode. Therapy is synchronized to a sensed ventricular event.

Progressive Episodes Therapies – This feature causes the device to skip therapies or modify high-voltage energy levels to ensure that each therapy delivered during an episode is at least as aggressive as the previous therapy.

Testing features

For more information about these features, see the reference manual.

Underlying Rhythm Test – This feature temporarily inhibits the pacing output of the device to enable the clinician to evaluate the patient's intrinsic heart rhythm. During the test, the device is temporarily programmed to a nonpacing mode.

Pacing Threshold Test – This feature allows the clinician to determine the patient's pacing stimulation thresholds. This information can be used to determine appropriate amplitude and pulse width settings that ensure capture and minimize output.

Wavelet Test – This feature evaluates the accuracy of the current Wavelet template and allows the clinician to collect a new template, if necessary.

Lead Impedance Test – This feature tests the integrity of the implanted lead system by measuring the impedance of the pacing and high-voltage electrodes. The test uses low-voltage, subthreshold pulses to make these measurements.

Sensing Test – This feature measures P-wave and R-wave amplitudes to help the clinician assess lead integrity and sensing performance. Mode, AV Delay, and Lower Rate can be programmed temporarily so that the device is not pacing the patient's heart, increasing the likelihood that sensed events will occur.

Charge/Dump Test – This feature tests the charge time of the capacitors and dumps any charge remaining on the capacitors.

EP Studies – This set of protocols allows clinicians to induce arrhythmias during electrophysiology studies. The available induction protocols are T-Shock, 50 Hz Burst, Fixed Burst, and Programmed Electrical Stimulation. Manual therapies are also available.

Additional operations

MRI SureScan feature – This feature allows patients to be safely scanned by an MRI machine when used according to the specified MRI conditions for use. Refer to the MRI technical manual for additional information.

Pacing mode information

Pacemaker modes are described using the NBG code. The five-letter NBG code, named after The North American Society of Pacing and Electrophysiology (NASPE) and the British Pacing and Electrophysiology Group (BPEG), describes the operation of implantable pulse generators. The NBG code, which supersedes the ICHD Code, is described in Table 1.

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Table 1. The Revised NASPE/BPEG Generic Code for antibradycardia pacing

Position:	I	II	III	IV	V
Category:	Chamber(s) paced	Chamber(s) sensed	Response to sensing	Rate modulation	Multisite pacing ^a
	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None A = Atrium V = Ventricle D = Dual (A + V)	O = None T = Triggered I = Inhibited D = Dual (T + I)	O = None R = Rate modulation	O = None A = Atrium V = Ventricle D = Dual (A + V)
Manufacturers' designation only:	S = Single ^b (A or V)	S = Single ^b (A or V)			

^a Medtronic device do not use the multisite pacing code.

^b The programmer displays A or V (not S) for chambers paced and sensed.

DDDR and DDD modes – When pacing in the DDDR or DDD mode, if the device senses an intrinsic atrial event, it schedules a ventricular paced event in response. The delay between the sensed atrial event and the corresponding ventricular paced event is the programmed Sens ed AV interval. The delay between the paced atrial event and the corresponding ventricular paced event is the Paced AV interval. If the current pacing interval ends before the device senses an atrial event, the device paces the atrium and then schedules a ventricular paced event after the Paced AV interval. If a ventricular sensed event occurs during the Sens ed AV interval or the Paced AV interval, ventricular pacing is inhibited. DDDR pacing occurs at the sensor rate, unless the intrinsic atrial rate is faster. DDD pacing occurs at the programmed lower rate, unless the intrinsic atrial rate is faster.

DDIR and DDI modes – When the device is pacing in the DDIR or DDI mode, sensed atrial events are not tracked. If the device senses an intrinsic atrial event, atrial pacing is inhibited, but a Sens ed AV interval is not started. Instead, ventricular pacing is delivered at the current pacing rate. If the current pacing interval ends before the device senses an atrial event, the device paces the atrium and then schedules a ventricular paced event after the Paced AV interval. If a ventricular sensed event occurs during the Paced AV interval, ventricular pacing is inhibited. DDIR pacing occurs at the sensor rate. DDI pacing occurs at the programmed lower rate.

DOO mode – The DOO mode provides AV sequential pacing at the programmed lower rate with no inhibition by intrinsic events. In the DOO mode, no sensing occurs in either chamber.

VVIR and VVI modes – In the VVIR and VVI modes, the ventricle is paced if no intrinsic ventricular events are sensed before the current pacing interval ends. Pacing occurs at the sensor rate in the VVIR mode and at the programmed lower rate in the VVI mode.

VOO mode – The VOO mode provides ventricular pacing at the programmed lower rate with no inhibition by intrinsic ventricular events. In the VOO mode, no ventricular sensing occurs.

AAIR and AAI modes – In the AAIR and AAI modes, the atrium is paced if no intrinsic atrial events are sensed before the current pacing interval ends. Pacing occurs at the sensor rate in the AAIR mode and at the programmed lower rate in the AAI mode.

AOO mode – The AOO mode provides atrial pacing at the programmed lower rate with no inhibition by intrinsic atrial events. In the AOO mode, no atrial sensing occurs.

MVP mode – MVP (Managed Ventricular Pacing) is an atrial-based pacing mode that is designed to switch to a dual chamber pacing mode during AV block.

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

	MVP allows the device to function in the AAIR or AAI pacing mode when AV conduction is intact but to switch to DDDR or DDD pacing mode during AV block. It also allows the device to perform periodic conduction checks, with the ability to switch back to AAIR or AAI mode when AV conduction resumes. MVP provides backup ventricular support for transient loss of AV conduction. The AAIR<=>DDDR mode switches between the AAIR and DDDR modes and the AAI<=>DDD mode switches between the AAI and DDD modes. ODO mode – The ODO mode does not deliver ventricular or atrial pacing outputs, regardless of the intrinsic rate. The ODO mode is intended only for those situations where the clinician wants to turn off bradycardia pacing outputs from the device.													
4. Conditions for use	In medical and preventive treatment facility. The cardioverter defibrillators should be used only by or under the supervision of physicians trained.													
5. Composition and list of parts and accessories	Digital dual chamber implantable cardioverter defibrillator Evera MRI S DR SureScan, model DDMC3D1: -1 implantable cardioverter defibrillator -1 torque wrench -1 DF-1 pin plug -Product Documentation* - CD-ROM *Product documentation includes - Device Manual, Limited warranty, Medical Procedure and EMI Precautions, Explanation of Symbols for ICD and CRT-D devices, Document: Electronic instructions for use/manuals, Appendix to the Device Manual (not in package content, printed and supplied locally on the territory of Russian Federation to the customers), Stickers, Reference manual (not in package content, available on http://manuals.medtronic.com/)													
6. Technical details of medical product	Overview of Physical Characteristics of the Evera MRI models: <table><tr><td>Evera MRI Surescan S DR</td><td>DDMC3D1</td></tr></table> Physical characteristics <table><tr><td>Volume^a</td><td>33 cm³ ± 0.5 cm³</td></tr><tr><td>Mass</td><td>77 g ± 1 g</td></tr><tr><td>H x W x D</td><td>66 mm (±1.0 mm) x 51 mm (±1.0 mm) x 13 mm (± 0.5 mm)</td></tr><tr><td>Surface area of device can</td><td>57 cm² ± 0.5 cm²</td></tr><tr><td>Radiopaque ID^b</td><td>PFZ*</td></tr></table>		Evera MRI Surescan S DR	DDMC3D1	Volume ^a	33 cm ³ ± 0.5 cm ³	Mass	77 g ± 1 g	H x W x D	66 mm (±1.0 mm) x 51 mm (±1.0 mm) x 13 mm (± 0.5 mm)	Surface area of device can	57 cm ² ± 0.5 cm ²	Radiopaque ID ^b	PFZ*
Evera MRI Surescan S DR	DDMC3D1													
Volume ^a	33 cm ³ ± 0.5 cm ³													
Mass	77 g ± 1 g													
H x W x D	66 mm (±1.0 mm) x 51 mm (±1.0 mm) x 13 mm (± 0.5 mm)													
Surface area of device can	57 cm ² ± 0.5 cm ²													
Radiopaque ID ^b	PFZ*													

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Medtronic radiopaque identifier^b

Materials in contact with human tissue^c

Titanium, polyurethane, silicone rubber

Battery

Hybrid CFx lithium/silver vanadium oxide


^a Volume with connector ports unplugged.

^b The radiopaque ID, which includes a Medtronic-identifier symbol, can be viewed in a fluoroscopic image of the device.

^c These materials have been successfully tested for the ability to avoid biological incompatibility. The device does not produce an injurious temperature in the surrounding tissue during normal operation.

* Radiopaque ID (PFZ) is made from tungsten, therefore it does not require testing. Tungsten is heavy metal and is visible on any X-rays.

Connector ports and suture holes


- 1 DF-1 connector port, SVC
- 2 DF-1 connector port, RV
- 3 Device Active Can electrode
- 4 IS-1 connector port, RV


Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

5 IS-1 connector port, A
6 Suture holes

Device radiopaque symbols: DDMC3D1

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Electrical specifications

Basic battery characteristics and device specification:

Battery characteristics	
Manufacturer	Medtronic Energy and Component Center
Model	M945899A
Chemistry	Hybrid CFx lithium/silver vanadium oxide
Battery electrical specifications	
Nominal voltage	3.2 V
Mean capacity to RRT	1.0 Ah
Minimum capacity after RRT	0.1 Ah
Device electrical specifications	
Pacing rate limit (protective feature) ^a	200 min ⁻¹
Input impedance	150 kΩ minimum

^a Does not apply during ATP therapies or ventricular safety pacing.

Peak ICD output voltage during high-voltage shock delivery:


Programmed energy	Peak voltage for first pulse phase	Peak voltage for second pulse phase
Minimum: 0.4 J (±0.25 J)	76 V (±16%)	36 V (±30%)
Mean: 18 J (±20%)	497 V (±10%)	249 V (±25%)
Maximum: 35 J (±15%)	699 V (±10%)	351 V (±25%)

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Typical output waveform shapes

1. Pacing waveform;
2. Monophasic high-voltage waveform (T-Shock inductions only), 50% tilt;
3. Biphasic high-voltage waveform, 50% tilt.


Measuring methods

Device parameters, such as pulse duration, pulse amplitude, and sensitivity, are measured at the standard conditions of $37^{\circ}\text{C} \pm 2^{\circ}\text{C}$ and $500\ \Omega \pm 1\%$ load according to EN 45502-2-1.

Pulse duration – Pulse duration is measured at 1/3 peak voltage levels according to standard EN 45502-2-1 (see Figure below). When applying this measurement method, the measured pulse width W depends on the load R_{load} (in ohms) and programmed pulse width W_p (in seconds) with tolerance $W \leq W_p + 34\ \mu\text{s}$ and $W \geq$ the smaller of $(W_p - 16\ \mu\text{s})$ or $[124\ \mu\text{s} + (4\ \mu\text{s} \times R_{\text{load}})]$.

Amplitude – The pulse amplitude is calculated according to standard EN 45502-2-1 (see Figure below). When applying this measurement method, the measured amplitude A depends on the programmed amplitude A_p and programmed pulse width W_p : $A = A_p \times [0.9 - (W_p \times 0.145\ \text{ms}^{-1})]$. The tolerance (+40%/-30% for voltages less than 2.0 V, and $\pm 30\%$ for voltages greater than or equal to 2.0 V) is applied not to the programmed setting, but to the calculated amplitude A .

Figure. Measurement of pulse duration and amplitude


1. Maximum amplitude;
2. 1/3 maximum amplitude;
3. Pulse duration (D);
4. Time integral over voltage (F);
5. Pulse amplitude (A).

Sensitivity – Sensitivity is defined as the voltage amplitude of a test signal that is just sufficient to be sensed by the device. The programmable values for RV sensitivity assume a 40 ms sine² waveform. When using the test signal defined in EN 45502-2-1 (see Figure below), the rated ventricular sensing threshold will be 1.5 times the programmed RV Sensitivity value.

Figure. Measurement of sensitivity

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Amplitude

Notes:

- When measuring the pacing and sensing parameters with pacing system analyzers, considerable differences may be observed with the specifications presented in this manual, because the measuring methods employed by such systems may differ from those described previously.
- Lead impedance measurement results may be distorted by electrocardiogram monitoring equipment.

Common mode rejection ratio – The common mode rejection ratio (CMRR) for frequencies 16.6 Hz, 50 Hz, and 60 Hz is at least 100 (40 dB). The calculation of the CMRR ratio was done based on measurements performed with the sinusoidal waveform injected directly into the device. The device and lead system CMRR ratio depend on several factors, such as the placement of the electrodes or electrode separation, and may be lower than the device CMRR ratio.

Variation with temperature

Basic rate, test pulse rate, pulse duration, and pulse amplitude remain within expected tolerances when the device temperature is between 22 °C and 45 °C. Sensitivity at nominal conditions as measured at 37 °C can vary as much as $\pm 1\%$ per °C from 22 °C to 45 °C.

Replacement indicators

The Remaining Longevity estimate, replacement status, and battery voltage appear on the programmer display and on printed reports. The Recommended Replacement Time (RRT) and the End of Service (EOS) conditions are listed in Table.

Replacement indicators

Recommended Replacement Time (RRT)	<2.73 V on 3 consecutive daily automatic measurements
End of Service (EOS)	3 months after RRT

Remaining Longevity – The Remaining Longevity estimate displays the estimated time remaining until device RRT.

RRT (Recommended Replacement Time) – The programmer displays the RRT battery status to indicate that replacement of the device is recommended.

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

RRT date – The programmer displays the date when the battery reached RRT on the Quick Look II and Battery and Lead Measurements screens.

EOS (End of Service) – The programmer displays the EOS battery status to indicate that the device should be replaced immediately and may not operate per specifications.

Replace at EOS – If the programmer indicates that the device is at EOS, replace the device immediately.

Prolonged Service Period – The Prolonged Service Period (PSP) is the time between the RRT and EOS. The PSP is defined as 3 months assuming the following conditions: 100% VVI pacing at 60 min⁻¹, 2.5 V RV pacing amplitude; 0.4 ms pulse width; 600 Ω pacing load; and 6 full-energy charges. The EOS may be indicated before the end of 3 months if the device exceeds these conditions.

Projected service life

The projected service life in years for the device is shown in Table 9. The data is based on pacing outputs programmed to VVI mode, specified pacing amplitude, 0.4 ms pulse width, and 60 min⁻¹ pacing rate.

The service life projections are based on the following assumptions:

- Semi-annual maximum energy charging frequency
- Pre-arrhythmia EGM storage programmed to On for a 6-month period (two 3-month follow-up intervals), over the entire life of the device
- 3 hours of wireless telemetry during implant
- A quarterly schedule of Medtronic patient monitors remote transmissions
- 1 hour of in-office wireless telemetry annually
- Typical shelf storage time before implant

Projected service life estimates are based on accelerated battery discharge data and device modeling as specified. Do not interpret these values as precise numbers.

Projected service life in years

Pacing mode, percent pacing	Pacing amplitude	Projected service life in years	
		Pacing impedance 500 Ω	Pacing impedance 600 Ω
DDD, 0%	2.5 V	9.7	9.7
	3.5 V	9.6	9.6
DDD, 15%	2.5 V	9.1	9.2
	3.5 V	8.7	8.8
DDD, 50%	2.5 V	8.3	8.5
	3.5 V	7.1	7.4
DDD, 100%	2.5 V	7.3	7.6
	3.5 V	5.7	6.1

Projected service life in years, with the Pacing Mode programmed to AAI<=>DDD

Projected service life in years

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Pacing mode, percent pacing	Pacing amplitude	Pacing impedance 500 Ω	Pacing impedance 600 Ω
AAI $\leftrightarrow$ DDD (MVP mode) 50% atrial/5% ventricular	2.5 V	8.9	9.0
	3.5 V	8.2	8.3

Projected service life considerations

Additional full-energy charges – Each additional full-energy charge due to therapy shock or device testing reduces projected service life by approximately 33 days.

Pre-arrhythmia EGM storage – Full-time use of Pre-arrhythmia EGM storage reduces projected service life by approximately 2.4 months per year, or 20%.

Medtronic patient monitor remote transmissions – Additional Medtronic patient monitor remote transmissions reduce projected service life. For example, a DR patient who received 0% atrial and ventricular pacing, DDD mode, with 600 Ω impedance would expect 9.7 years projected service life. Projected service life reductions for more frequent remote transmission rates are as follows:

- Monthly transmissions over the life of the device reduce projected service life by 80 days, or 2%.
- Weekly transmissions over the life of the device reduce projected service life by 381 days, or 11%.
- Daily transmissions over the life of the device reduce projected service life by 1549 days, or 44%.
- A single additional transmission reduces projected service life by approximately 0.8 days, or 0.02%.

Wireless telemetry – Each additional hour of wireless telemetry use (in-office or implant) reduces the projected service life by approximately 7.6 days, or 0.25%.

Shelf storage time – Maximum shelf storage time of 18 months reduces projected service life by approximately 6.4%.

Energy levels and typical charge times

Energy levels – Stored energy is always greater than the delivered energy. Stored energy is derived from the peak capacitor charge.

Typical charge times – The most recent capacitor charge time appears on the programmer display and on printed reports. You can evaluate charge time using the Charge/Dump Test.

Maximum energy levels and typical full energy charge times

Maximum programmed energy	35 J
Maximum delivered energy ^{a,b}	36 J
Maximum stored energy ^c	42 J
Typical charge time at Beginning of Service (BOS) ^d	8.3 s
Typical charge time at Recommended Replacement Time (RRT) ^d	12.3 s

^a Energy delivered at connector block into a 50 Ω load.

^b For 35 J programmed energy, delivered energy exceeds 35 J.

^c Energy stored at charge end on capacitor.

^d Charge time during a nonwireless telemetry session may be slightly higher.

Magnet application

When a magnet is placed near the device, tachyarrhythmia detection is suspended and no tachyarrhythmia therapies are delivered.

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Alert tones sound if programmed. The device ignores the magnet in the programmer head when telemetry communication is established through the programmer head. Before implant and for the first 6 hours after implant, the device does not sound audible tones when a magnet is placed over the device.


Note: If MRI SureScan is programmed to On, tachyarrhythmia detection and Medtronic CareAlerts (including audible alerts) are suspended.

Torque wrench (included in the package with ICD) – Use only the torque wrench supplied with the device. The torque wrench is designed to prevent damage to the device from overtightening a setscrew. Other torque wrenches (for example, a blue-handled or right-angled hex wrench) have torque capabilities greater than the lead connector can tolerate.

Materials (not in contact with the human body):

Thrust Bearing	Polytetrafluoroethylene per AMS 3652
Plastic Components	Polyetherimide (G.E. Ultem 1100F-8104); Color – Off-white, or equivalent
Hex Shaft	Carpenter Custom 455 Corrosion resistant steel – vacuum heat treated (passivated per Medtronic Specification – 200034-001)
Other Metal Components	Corrosion resistant steel (passivated per Medtronic Specification – 200034-001)

Dimensional (inches): (tolerance: .xx ± .02; .xxx ± .005)

The torque is limited to 12 to 17 in-oz (0.08 – 0.12 N*m). during first clockwise revolution; it will click when torque value in this range is achieved.


Hardness: 47-55 Rc

Mass: 8 g ± 0.5 g

DF-1 Connector Pin Plug

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential


DF-1 Connector Pin Plug (unit: mm)

DF-1 Connector Pin Plug included in package with the following models only: Evera MRI (DDMC3D1). This device is a plug to fill an IS-1 (ISO 5841-3) compliant connector port of pulse generators, such as an implantable cardiac pacemaker, to which an implantable IPG/ICD lead connector is not connected.

Materials: Handle – Silicone rubber (Silastic Biomedical Grade Liquid Silicone Rubber (LSR); Shaft – Stainless steel 316L

Mass: 0.41 gramm +/- 0.005 gramm

Insertion / With drawal – meets all ISO requirements 11318 (DF-1); ISO 5843 (IS-1)

The strength of the connection between the medical device and the electrode (introduction / extraction):

The connector socket accepts the corresponding pass gauge. The insertion / removal force of the connection socket must not exceed 2.0 lbs. Force (9 N) for DF-1 / IS-1 connectors.

IS-1, DF-1 Low Voltage Leakage Impedance - IS-1 and DF-1: Port leakage impedance shall be 50 kΩ minimum between conductive elements intended to be electrically insulated by sealing rings.

Mechanical test

The intent of mechanical vibration, mechanical shock, and controlled release free fall drop conditioning is to assure that the device performance is not adversely affected by the mechanical environment seen during shipping, distribution, and clinical pre-implant life of the device.

The 11-pin FT left hand shield shows the same stiffness (if not slightly lower – worse case) as the 9-pin shield according to Finite Element Analysis. The analysis compares deflections for the two left hand shield configurations. In the analysis, the shield is fixed around the entire edge and the load is applied similar to shield deflection test load, at the center of the shield. The output is the

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

deflection response for different force magnitudes and stress comparison at the FT shelf area. The distribution of stress at the FT shelf area is shown in see Figure below. Dark blue indicates a minimum value of stress. As stress increases, the color changes to light blue and then green, so the 11-pin FT shelf shows a higher stress distribution.



Figure 10: Stress Distribution on the Left Hand Shield (Left: 9-pin FT; Right: 11-pin FT)

The shield bracket location does not affect mechanical DAU testing of the Electronic Module Assembly (EMA). The connector attach is evaluated as part of the connector DAU testing and not the Mechanical DAU testing. Mechanical DAU is not designed to evaluate ICRs and Connector Attach.

Leak Rate

Each device is tested in manufacturing to assure an actual leak rate less than or equal to 1×10^{-8} standard cm³/sec of helium.

Application Software (External) for Evera

Safety class SW – C

The Model SW033 software is version 8.2

The model SW033 application software is loaded on the approved Medtronic Analyzer programmer and is responsible for providing support for the features and functionality for the device. The SW033 software includes all the executable and non-executable files and data files needed to support programming the Evera devices on the Medtronic Analyzer programmer systems. The software allows the clinician to interrogate and program the device, navigate through the user interface to interact with the implantable device, run tests, set up therapies, perform data analysis, and print reports.

7. The possibility and ways of integrating with other medical products

Leads – The lead system used with this device must provide sensing, pacing, and cardioversion and defibrillation therapies to the right ventricle (RV) and sensing and pacing to the atrium (A). Do not use any lead with this device without first verifying lead and connector compatibility.

The device typically is implanted with the following leads:

- 1 quadripolar/tripolar transvenous lead with a DF4-LLHH or DF4-LLHO connector in the right ventricle (RV) for sensing and pacing, and for cardioversion and defibrillation therapies
- 1 bipolar transvenous lead with an IS-1 connector in the atrium (A) for sensing and pacing. Use of a bipolar atrial lead with ring and tip electrodes spaced ≤ 10 mm apart to reduce far-field R-wave sensing is recommended. *Example of compatible leads produced by Medtronic Inc. – Sprint Quattro family*

Programmers and software – The Medtronic CareLink Programmer and software are used to program this device. Refer to the reference manual for information about using the programmer. Programmers from other manufacturers are not compatible with Medtronic devices, but they do not damage Medtronic devices.

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

	Model 27901 Conexus Activator – The Medtronic Model 27901 Conexus Activator allows you to turn on Conexus wireless telemetry for implanted devices that support wireless telemetry. The Conexus Activator is used with the Medtronic CareLink Model 2090 Programmer with Conexus telemetry in the hospital or clinic. Medtronic Analyzer – The system supports the use of the Medtronic Analyzer, an accessory of the Medtronic CareLink Programmer. The system allows you to have a device session and an analyzer session running at the same time, to switch quickly from one session to the other without having to end or restart sessions, and to send data from the analyzer to the programmer. Medtronic Patient Monitor – Patients use the Medtronic patient monitor to gather information automatically from their implanted devices and communicate the information to their physicians through the Medtronic CareLink Network. For information on using the patient monitor, refer to the patient manual; for connection and usage information, refer to the monitor literature. CareLink Network – The Medtronic CareLink Network is an Internet-based service that allows a patient to transmit full cardiac device information from home or other locations to the physician over a secure server. The CareLink Network may be unavailable in some geographic locations. CareLink Express monitor Model 2020B – The Medtronic CareLink Express monitor Model 2020B is a portable interrogation and data transfer tool used with Medtronic implantable cardiac devices. Model 2696 InCheck Patient Assistant – Patients can use the Medtronic Model 2696 InCheck Patient Assistant to initiate recording of cardiac event data in the device memory, verify whether the implanted device has detected a suspected atrial tachyarrhythmia, and, if the device is programmed to allow patient-activated cardioversion, request delivery of atrial cardioversion therapy. Model 2490C CareLink Monitor – Patients use the Medtronic Model 2490C CareLink Monitor to gather information automatically from their implanted devices and communicate the information to their physicians. For information on using the CareLink Monitor, refer to the patient manual; for connection and usage information, refer to the monitor literature.																		
8. EMC information	<table><tr><th colspan="3">Electromagnetic emission</th></tr><tr><td colspan="3">Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment</td></tr><tr><td>Electromagnetic emission test</td><td>Conformity</td><td>Electromagnetic environment – instructions</td></tr><tr><td>Radio interference according to CISPR 11</td><td>Group 1</td><td rowspan="4">Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.</td></tr><tr><td>Radio interference to CISPR 11</td><td>Class B</td></tr><tr><td>Harmonic components of the consumed current in accordance with IEC 61000-3-2</td><td>Do not use</td></tr><tr><td>Voltage fluctuations and flicker according to IEC 61000-3-3</td><td>Do not use</td></tr></table>	Electromagnetic emission			Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment			Electromagnetic emission test	Conformity	Electromagnetic environment – instructions	Radio interference according to CISPR 11	Group 1	Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.	Radio interference to CISPR 11	Class B	Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use	Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use
Electromagnetic emission																			
Digital Implantable dual chamber Cardioverter Defibrillators Evera are intended for use in the electromagnetic environment defined below. Buyer or user Digital Implantable dual chamber Cardioverter Defibrillators Evera should ensure its use in the specified environment																			
Electromagnetic emission test	Conformity	Electromagnetic environment – instructions																	
Radio interference according to CISPR 11	Group 1	Digital Implantable dual chamber Cardioverter Defibrillators Evera should emit electromagnetic energy to perform the main function. Digital Implantable dual chamber Cardioverter Defibrillators Evera are suitable for use in all locations.																	
Radio interference to CISPR 11	Class B																		
Harmonic components of the consumed current in accordance with IEC 61000-3-2	Do not use																		
Voltage fluctuations and flicker according to IEC 61000-3-3	Do not use																		

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Noise immunity			
Immunity test	The test level of IEC 60601	Level of compliance	Electromagnetic environment - instructions
Electrostatic Discharge (ESD) in accordance with IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%
Magnetic field of industrial frequency according to IEC 61000-4-8	3 A/m	0,3 A/m	The magnetic fields of the industrial frequency must be measured at the designated installation site to ensure that the field strength is sufficiently low
Immunity test	The test level of IEC 60601	The actual level of noise immunity	Level of compliance
Conductive noise induced by radio-frequency electromagnetic fields in accordance with IEC 61000-4-6	3 V (rms) in the band from 150 kHz to 80 MHz outside the frequencies allocated for the LPM of HF devices	1.7 V (rms value)	1 V (rms value)
The emitted radio-frequency electromagnetic field according to IEC 61000-4-3	10 V (RMS) in the band from 150 kHz to 80 MHz at the frequencies allocated for the LPM of high-frequency devices	2.3 V (rms value)	1 V (rms value)
		13 V/m	10 V/m
Noise immunity			
The test for Noise immunity	Test Level according to IEC 60601	Level Conformity	Electromagnetic environment - instructions
Electrostatic Discharge (ESR) according to IEC 61000-4-2	± 6 kV - contact discharge ± 8 kV - air discharge	± 6 kV - contact discharge ± 8 kV - air discharge	The floors of the room must be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30%.

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	A magnetic field Industrial Frequency (50/60 Hz) IEC 61000-4-8	3 A/m	3 A/m	The levels of the magnetic field of the industrial frequency must correspond to the typical conditions of a commercial or hospital environment.
	Radio frequency Electromagnetic field IEC 61000-4-3	3 V / m In the band from 80 MHz to 2.5 GHz	3 V/m	However, the readings of the weights built into the bed do not have a critical value for the operation of a medical device: «Digital Implantable dual chamber Cardioverter Defibrillators Evera» $D = 1,2\sqrt{P}$ in the band from 26 to 800 MHz $D = 2,3\sqrt{P}$ in the band from 800 MHz to 2.5 GHz, where P is the nominal maximum output transmitter power in Watts (W) set by the manufacturer, and d - Recommended spatial clearance in meters (m). The field strength in the propagation of radio waves from stationary radio transmitters from the results of observations of the electromagnetic situation should be lower than the level of compliance in each frequency band. Interference may occur near equipment  marked with a sign. This device may not cause harmful interference, and this. To accept any interference received, including interference that may cause undesired operation of the device. This device must not cause harmful interference to radio stations operating in the 400.150-406.000 MHz band and used for meteorological purposes, for meteorological satellites and Systems of Earth exploration satellites. This device must accept any interference received, including interference that may cause undesired operation.

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

NOTE 1 Higher frequency bands are applied at 80 MHz and 800 MHz.

NOTE 2: This manual is not applicable in all situations. The propagation of electromagnetic waves is affected by absorption and reflection from structures, objects and people.

9. Biocompatibility Testing

Table of Tissue Contacting Materials:

Material	Device Components
Titanium grade 1 <u>ASTM F67</u> Fe to 0.12 C to 0.05 Si to 0.08 N to 0.03 Ti 99.5 – 99.99 O to 0.1 H to 0.008	LH shield
	RH shield
	Connect attach pin
	Fastener bracket
Titanium grade 2 <u>ASTM F67</u> Fe to 0.18 C to 0.07 Si to 0.1 N to 0.03 Ti 98.61 – 99.7 O to 0.12 H to 0.01	End cap
Silicone rubber MED 2000	Adhesive
Polyurethane 75D Elasthane 75D™/Tecothane TT-1075D-M	Connector module
Silicone rubber ETR Elastomer Q7-4735	Split grommet
LSR with titanium dioxide 2% (Q74850, colorant Ferro white 61-182000)	IS4 seal

Evera (Blackwell) has been evaluated as guided by the applicable sections of ISO 10993-1:2009 that pertain to biological effects. Table above lists all the subparts of the ISO 10993 standard and how they relate to this evaluation. The materials used in Evera (Blackwell) have a long history of safe clinical use in previous Medtronic devices with equivalent patient tissue contact applications and manufacturing processing, and extensive available biological safety test data. ISO 10993-1:2009 section 4, General principles applying to biological evaluation of medical devices, states: "Evaluation may include both a study of relevant preclinical and clinical experience and actual testing. Such an evaluation may result in the conclusion that no testing is needed if the material has a demonstrable safe history of use in a specified role and physical form that is equivalent to that of the device under design."

Medtronic has performed extensive biological safety testing on each material separately and has demonstrated a safe history of

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

clinical use of each material in predicate market-approved devices. The testing was performed in the past on other Medtronic parts made of the same materials. The biological safety test data and safe history of use are summarized in each material's appendix attached to this evaluation. In instances where the testing was performed according to previous versions of standards, gap analyses between the standards used and the current ISO 10993 standards are documented in the appendices to demonstrate the applicability of the testing. Based on Table 3.1 below, the documented safe history of clinical use, the biological safety testing performed, and the gap analyses between the standards, all tissue-contacting materials used in Evera (Blackwell) are compliant with the current revisions of all applicable ISO 10993 series standards.


In the case of the new 75D polyurethane connector modules, this component undergoes additional processing compared to the same material used in predicate devices. Therefore, focused chemical characterization and biological safety testing was performed on the new Evera (Blackwell) connector module.

The history of use of the materials in predicate devices, the existing biological safety test data available for the materials, and additional biological safety testing, chemical characterization, and toxicological risk assessment showing no additional biological safety concerns resulting from the connector module processing, demonstrate acceptable biological safety risk for the Evera (Blackwell) device as guided by ISO 10993-1:2009.

A biological evaluation as guided by the applicable sections of ISO 10993-1:2009 that pertain to biological effects has been performed for Evera (Blackwell). Existing information demonstrates acceptable biological safety risk for the materials used in Evera (Blackwell). Additional biological safety testing, chemical characterization, and toxicological risk assessment was conducted as outlined in the biological evaluation plan.

10. Design and appearance, photo or image, the General Scheme

Evera S DR MRI DDMC3D1


Torque wrench and DF-1 pin plug (for Evera S DR MRI DDMC3D1):

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Product documentation Evera S DR MRI DDMC3D1

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Not to be reproduced without written permission from Medtronic, Inc.
 Medtronic, Inc.
 - Confidential

11. Information on labelling, packaging	Products supplied in packaging marked with the following information: • product name and catalog number; • LOT of the product; • trademark and address of the manufacturer; • single-use; • sterile, do not resterilize; • label complies with European Council 90/385/EEC; • manufacturing date and shelf life of the product; • storage conditions; • do not use if the package is damaged; • characteristics of the product
	Evera S DR MRI DDMC3D1

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential



Medtronic

Evera MRI™ S DR SureScan™ DDMC3D1

DDE-DDDR

MR Conditional

Defibrillation

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Burst RV

Burst RV

Burst RV

Burst RV

Burst RV

Burst RV

Adaptive

AAI ↔ DDD

↑ 130 min⁻¹
↓ 60 min⁻¹
Upper tracking rate and lower rate

0.5 mV
Sensitivity I/A

0.3 mV
Sensitivity RV

→ 180 ms
→ 160 ms
AV delay (post-ventricular)

± 3.5 V ±
→ 0.4 ms ±
Amplitude and pulse width NA

± 3.5 V ±
→ 0.4 ms ±
Amplitude and pulse width RV

OFF
Detects VT, FVT, VF

OFF
Therapy (RV, VT, FVT)

35/42 J
Threshold (sensed and stored) V

ON
Active Can™

OFF
Inhibits A1, A2

OFF
Inhibits A1, A2

35°C
131°F

PDO BY LIFETIME IMPLANTS
IC: 34080 IMPLANTS

STERILE EO

SN WWW000000W

2000-00-00
YYYY-MM-DD

2000-00-00
YYYY-MM-DD

www.medtronic.com/mimipulse

Medtronic, Inc., 710 Medtronic Parkway
Minneapolis, MN 55432 USA

Manufactured in
Tetechenat, Switzerland
www.medtronic.com

CE0123


DR00000001 A - 000000001 A

Sterile Pack Label

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

DF-1 Evera MRI™ S DR SureScan™ DDMC3D1 DDE-DDDR

Active

AAI <=> DDD

↑ 130 min⁻¹
↓ 80 min⁻¹

Upper tracking rate and lower rate

0.3 mV
Sensitivity: RA

0.3 mV
Sensitivity: RV

→ 180 ms
→ 150 ms

A-V interval (paced and sensed)

1.3 V E_T
→ 0.4 ms E_T
Amplitude and pulse width: RA

1.3 V E_T
→ 0.4 ms E_T
Amplitude and pulse width: RV

OFF
Detection: VI, FVI, VF

OFF
Threshold: (RV), VI, FVI

55442 J
Threshold (diagnosed and treated): VF

ON
Active Can™

OFF
Threshold: AT, AF

OFF
Threshold: AT, AF

05 °C
131 °F
-10 °C
0 °F
STERILE ED
DO NOT USE FOR IMPLANTS
DO NOT IMPLANT

Medtronic Inc., 711 Medtronic Parkway
Minneapolis, MN 55424 USA
Medtronic B.V., Earl Beelaersstraat 10
6422 PZ Venlo, The Netherlands
www.medtronic.com

Manufactured in:
Tuesdun, Switzerland
Product of Switzerland

SN WWW000000W
2000-00-00
YYMMDD
2000-00-00
YYMMDD

CE0123
Medtronic A + B (Europe) A

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
• Confidential

Medtronic

2x 2x

Evera MRI™ S DR SureScan™ DDMC3D1

(01)00000000000000(17)000000(21)WWW000000W

REF DDMC3D1
SN WWW000000W
2000-00-00

MR5555A001 A - MR5512A001 A

Medtronic

2x 2x

- PhysioCurve™ Design
- SmartShock™ Technology

MR5555A001 A - MR5512A001 A

REF DDMC3D1
SN WWW000000W
2000-00-00
PIN: WWWWWW00W

Shelf Box Label

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Medtronic A001 A - Medtronic A001 A

2000-00-00
YYYY-MM-DD

Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ SN WWW000000W

Medtronic
Evera MRI™ S DR SureScan™ DDMC3D1
2000-00-00
YYYY-MM-DD
REF DDMC3D1
SN WWW000000W


Medtronic
Evera MRI™ S DR SureScan™ DDMC3D1
2000-00-00
YYYY-MM-DD
SN WWW000000W

Medtronic
Evera MRI™ S DR SureScan™ DDMC3D1
2000-00-00
YYYY-MM-DD
SN WWW000000W

FYR Labels

Sticker for final package will be added locally in Russia after approval with information below (information can be changed in accordance with Russian legislation):

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

 	Product name (in Russian)	
	Manufacturer (name)	UPN (unique product number)
	Country of origin	CFN n (customer facing number)
	See Device Manual	
	RC # _____ date of issue _____	
		Conformite Européenne (European Conformity). This symbol means that the device fully complies with AIMD Directive 90/385/EEC (NB 0123) and R&TTE Directive 1999/5/EC
		This symbol means that the device fully complies with the Australian Communications and Media Authority (ACMA) and the New Zealand Ministry of Economic Development Radio Spectrum Management standards for radio communications products
		Reorder number
		Serial number
		Manufacturer
		Sterilized using ethylene oxide

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

		Use by
		Do not use if package is damaged
		Date of manufacture
		Do not reuse
		Authorized representative in the European Community
	AAI ⇌ DDD	AAI ⇌ DDD mode
		Package contents
		Torque wrench
		Pin plug, DF-1
		Product documentation

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

		Implantable device
		ICD (dual chamber)
		Header shape
		DF-1 connector
		IS-1 connector
		Adaptive
		Upper tracking rate and lower rate
		Open here
		www.medtronic.com/manuals
		Temperature limitations

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

		Dangerous voltage
		Amplitude and pulse width: RA
		Amplitude and pulse width: RV
		Sensitivity: RV
		Sensitivity: RA
		A-V interval (paced and sensed)
		Therapies: AT, AF
		Intervention: AT, AF
		Detection: VT, FVT, VF
		Therapies (RV): VT, FVT

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

		Therapies (delivered and stored): VF
		Active Can
		Defibrillation
		AV cardioversion
		Ramp: RA, RV
		Ramp+: RV
		Burst: RV
		Burst+: RA
		50 Hz Burst
		SureScan symbol

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	MR Conditional (Evera MRI SURESCAN XT DR Evera MRI SURESCAN S DR)
	PRODUCTS and SYSTEMS that include radio frequency transmitters or who intentionally use radio frequency electromagnetic energy for diagnosis or treatment shall be marked in the form of a sign of non-ionizing radiation
	This conformity mark is placed on the products for which the mandatory certification procedure was carried out

Test performance - packaging in accordance with ISO 11607-1.

Environmental Conditioning

- The package must not allow the contents to spill out or be separated from the package.
- Glue joints, tape joints, tape and overwraps applied to the package must remain attached and if printed must remain legible.
- Labels must remain attached and legible.
- Label lift may occur around the edges of label stock, and bubbles in the interior of the label are expected. These conditions are acceptable as long as the label is still attached and legible.
- The package must remain closed and intact per design intent.

Shelf Box/Protective Package

- Boxes and cartons may be creased and bent and corners crushed, cracked or split; Shrink wrap may be wrinkled, torn and may have split seals.
- Cushioning may be cracked or split; compressed and dented.
- The shrink wrap must remain on the package and there must be evidence of a tear tape or microperforation.
- Seal labels or closure tape must keep the shelf box/protective package closed. Partial tears of seal labels and closure tape is acceptable.

Labeling

- Label lift may occur around the edges of label stock, and bubbles in the interior of the label are expected. These conditions are acceptable as long as the label is still attached and legible.
- Label stock may be torn or split.
- Removable adhesive labels must be removable without damaging the substrate or delaminating the label stock.

Sterile Barrier System

- Wrinkling of lids or flexible pouches.
- Any migration of the product or accessories within the package shall be reviewed and documented in the report. Displaced snap-in accessories inside blister packages such as wrenches or vein lifters are acceptable as long as they are prevented

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

from free circulation in the package which could cause device damage or compromise the sterile barrier.

- The seal flanges on a blister package may be bent.
- The sterile barrier package must also be able to be opened by the end user using good aseptic handling techniques.

All requirements were met during packaging tests

Technical characteristics of packaging: *

1. Microbial barrier:

- Maximum possible permeability of substrate: 5% Pmax

2. Toxicological properties:

Primary packaging:

- Sterile tray: Pre-shaped (thermoformable) polyethylene terephthalate glycol (PETG) copolymer without coating.

- Cover: tightly closed with a lid of fibrous Tyvek material on the basis of polyethylene, using the "spunbond" technology, with the adhesive composition Oliver 18B

Secondary packaging:

- Box: Polypropylene, molded under pressure

- Heat Shrink Film: Transparent polyalkene

3. Physical and chemical properties:

Overall dimensions: 16.9 cm $\pm$ 1 cm $\times$ 13 cm $\pm$ 0.5 cm $\times$ 2.1 cm $\pm$ 0.4 cm

Film thickness: 0.25 mm $\pm$ 0.05 mm

Density: 1.27 g / cm³

Transparency: 0.5%

Gloss level: 108 gloss units, at an angle of 45°

Spectral transmittance: 90%

Stress at the point of fluidity: 48 MPa

Young's modulus: 1720 MPa

Coefficient of friction: > 1

Coefficient of transmission of water vapor: 6 g / m²

3.1 Characteristics of Tyvek fibrous material:

Bendtsen air permeability: 572 ml / min

Permeability rate of water vapor: 1615 g / m²

Height of hydrostatic pressure: 1442 Pa

Tensile strength in the longitudinal direction: 196 N

Tensile strength in transverse direction: 200 N

Elongation in the longitudinal direction: 20%

Relative elongation in the transverse direction: 24%

Elmendorf fracture test in the longitudinal direction: 3.3 N

Elmendorf fracture test in the transverse direction: 3.5 N

Resistance to tearing along the Mullen: 1213 kPa

Impact strength according to Spencer: 8756 J / m²

Opacity: 91%

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Thickness: 178 microns

4. Compatibility with the methods of sterilization with which the use of the material is intended:

Methods of sterilization:

- Ethylene oxide
- Radiation sterilization

5. Compatibility with the processes of forming and gluing:

Termoskleivaetsya with a large number of thermoplastics, as PVC, styrene, PETG, HIPS; With polyethylene some kinds of polypropylene.

6. Shelf life after sterilization: 18 months

7. The materials are non-alkaline, non-toxic and odorless.

8. There are no holes, cracks, tears, wrinkles and clarifications in the materials.

9. Density: 1.27 g / cm³

10. The level of purity of the material: acceptable.

11. Specific and minimum physical characteristics:

Tensile strength (tear resistance): 1.25 N ± 0.1 N

Tensile load (tensile strength): 50 MPa

Thickness deviation: ± 0.05 mm

Breathability:

- for CO₂: 30 cm² / s-kgs / cm²
- for O₂: 10 cm² / s-kgf / cm²
- for N₂: 5 cm² / s-kgs / cm²
- Impact resistance (bursting strength):

At 23 ° C: 425 g ± 30 g

At -29 ° C: 350 g ± 25 g

12. pH: Tyvek has a neutral pH of 7.

13. Adhesive characteristics:

Peel strength (strength of the adhesive layer): 4.53 H-18.14 H

Weight of the adhesive layer: 4.5 kg ± 0.7 kg / 279 m²

Adhesion transfer: > 75%

Shelf life: <48 months

White color

The density of glue application: > 1.03 (g / cm³)

14. Components of packages do not react, do not contaminate, do not pass, do not affect negatively on products before, during or after sterilization.

15. Width of the adhesive layer: 10 mm ± 1 mm

16. Process indicators printed on packages meet the requirements of ISO 11140-1.

17. The opening of the packages by delamination is continuous and homogeneous, without delamination / rupture of the material.

18. Used dyes:

Inks are black C70, produced by Sun Chemical Corp., USA;


Blue SUNSPERSE 15, produced by Sun Chemical Corp, USA

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

Packaging Description
Shelf Box - Secondary Packaging

Package Part	Material
Box	Molded polypropylene

Not to be reproduced without written permission from Medtronic, Inc.
 Medtronic, Inc.
 - Confidential

Shrink Wrap


Clear olefin

Sterile package:

Package dimensions:


1. Sterile package lid (in inches)

Inner package lid:

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

Outer package lid:

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	<table> <tr> <th>Package Part</th><th>Material</th></tr> <tr> <td>Sterile tray</td><td>Pre-shaped (thermo-formed) uncoated PETG (polyethylene terephthalate, glycol modified)</td></tr> <tr> <td>Lid</td><td>Spun bonded olefin lids coated with Oliver 18B adhesive</td></tr> </table> * The deviation of the technical characteristics of the package is $\pm 5\%$, unless otherwise specified.	Package Part	Material	Sterile tray	Pre-shaped (thermo-formed) uncoated PETG (polyethylene terephthalate, glycol modified)	Lid	Spun bonded olefin lids coated with Oliver 18B adhesive														
Package Part	Material																				
Sterile tray	Pre-shaped (thermo-formed) uncoated PETG (polyethylene terephthalate, glycol modified)																				
Lid	Spun bonded olefin lids coated with Oliver 18B adhesive																				
12. Methods and conditions of sterilization time maintaining sterility (for sterile medical devices produced)	<table> <tr> <th colspan="2">Sterilization&Validation</th></tr> <tr> <td>Sterilization Method</td><td>Ethylene Oxide (EO)</td></tr> <tr> <td>Sterilization cycle/Process name</td><td>30-minute 100% EO sterilization process</td></tr> <tr> <td>Applicable Sterilization Process Validation Standard</td><td>ISO 11135-1:2007</td></tr> <tr> <td>Sterilant Gas type</td><td>100% EO</td></tr> <tr> <td>Relative Humidity (RH)</td><td>55.0-85.0%</td></tr> <tr> <td>EO Gas Concentration</td><td>≥ 720 mg/l</td></tr> <tr> <td>Temperature during EO Gas exposure</td><td>45.0-55.0 °C</td></tr> <tr> <td>Pressure during EO gas Exposure</td><td>55.0-75.0 kPa</td></tr> <tr> <td>EO Gas Exposure time</td><td>27-33 min</td></tr> </table> Medtronic has sterilized the package contents with ethylene oxide before shipment (30 min. 100 % EO). A SAL in excess of 10^{-6} is achieved. Detected EO residual levels were less than 4 mg and detected ECH residual levels were less than 9 mg. Therefore, compliance with ISO 10993-7.	Sterilization&Validation		Sterilization Method	Ethylene Oxide (EO)	Sterilization cycle/Process name	30-minute 100% EO sterilization process	Applicable Sterilization Process Validation Standard	ISO 11135-1:2007	Sterilant Gas type	100% EO	Relative Humidity (RH)	55.0-85.0%	EO Gas Concentration	≥ 720 mg/l	Temperature during EO Gas exposure	45.0-55.0 °C	Pressure during EO gas Exposure	55.0-75.0 kPa	EO Gas Exposure time	27-33 min
Sterilization&Validation																					
Sterilization Method	Ethylene Oxide (EO)																				
Sterilization cycle/Process name	30-minute 100% EO sterilization process																				
Applicable Sterilization Process Validation Standard	ISO 11135-1:2007																				
Sterilant Gas type	100% EO																				
Relative Humidity (RH)	55.0-85.0%																				
EO Gas Concentration	≥ 720 mg/l																				
Temperature during EO Gas exposure	45.0-55.0 °C																				
Pressure during EO gas Exposure	55.0-75.0 kPa																				
EO Gas Exposure time	27-33 min																				
13. Requirements for the protection of the environment in the application of medical products	Follow local regulations for proper disposal of this product. Do not dispose of this product in the unsorted municipal waste stream. This product contains materials that can harm the environment. See http://recycling.medtronic.com for instructions on proper disposal of this product. Explant and disposal Consider the following information related to device explant and disposal: • To prevent the device from delivering unwanted shocks, interrogate the device and disable tachyarrhythmia detection before explanting, cleaning, or shipping the device. • Explant the implantable device postmortem. In some countries, explanting battery-operated implantable devices is mandatory because of environmental concerns; check the local regulations. In addition, the device may explode if subjected to incineration or cremation temperatures. • Medtronic implantable devices are intended for single use only. Do not resterilize and reimplant explanted devices. • Contact Medtronic for Return Mailer Kits to return explanted devices for analysis and disposal. See the back cover for addresses.																				

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	Note: Disposal of explanted devices or leads is subject to local, state, and federal regulations. Waste class – B (in accordance with Russian classification).
14. Information about manufacturer of medical product (name, location address, mailing address, phone number, address of place of production)	Design facility: Medtronic Inc., USA; 8200 Coral Sea Street Ne, Mounds View, MN 55112, USA; Tel.: +1 763 526 2388; www.medtronic.com Manufacturer: Medtronic Inc., USA; 710 Medtronic Parkway N. E. Minneapolis, MN 55432, USA; Tel.: +1 763 514 4000, Fax: +1 763 514 4879; www.medtronic.com Representative of Manufacturer: «Medtronic» LLC; 123317, Moscow, Presnenskaya Naberezhnaya, 10; Tel.: (495) 580-73-77, Fax: (495) 580-73-78; www.medtronic.ru Legal person in whose name may be issued the certificate of registration (name, address, tel, fax, e-mail, website): «Medtronic» LLC; 123317, Moscow, Presnenskaya Naberezhnaya, 10; Tel.: (495) 580-73-77, Fax: (495) 580-73-78; www.medtronic.ru Addresses of manufacturing facility of medical device: 1. Medtronic Europe Sarl, Route du Molliau 31, Case Postale, 1131 Tolochenaz, Switzerland; 2. Medtronic Puerto Rico Operations Co., Juncos Road 31, Km. 24, Hm 4, Ceiba Norte Industrial Park, Juncos, PR 00777 USA.
15. Classification of medical product	Hazard class, depending on the possible risk associated with the use of the medical device - CLASS 3 / III; GMDN code – 37265 Dual-chamber implantable defibrillator; Class of electrical safety, laser safety and other hazards (if available): not applicable to ICDs/IPGs; The software security class: C (according to IEC 62304); Type of contact with the body – Prolonged contact with blood, tissues and body fluids; Type of working part: similar to CF; Sterility: sterile; Classification, depending on the possible consequences of failure during the use of the product: A; Class of medical waste: B.
16. Methods and apparatus for	Sterile. For single use only. Products cannot be resterilized.

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

disinfection and sterilization cleaning																																			
17. Shelf life period	18 months Projected service life The projected service life in years for the device is shown in the tables below. The data is based on pacing outputs programmed to the specified mode, the specified pacing amplitude, 0.4 ms pulse width, and 60 min⁻¹ pacing rate. The service life projections are based on the following assumptions: • Semi-annual maximum energy charging frequency • Pre-arrhythmia EGM storage programmed to On for a 6-month period (two 3-month follow-up intervals), over the entire life of the device • 3 hours of wireless telemetry during implant • A quarterly schedule of Medtronic CareLink Monitor remote transmissions • 1 hour of in-office wireless telemetry annually • Typical shelf storage time before implant Projected service life estimates are based on accelerated battery discharge data and device modeling as specified. Do not interpret these values as precise numbers. Projected service life in years <table><tr><th rowspan="2">Pacing mode, percent pacing</th><th rowspan="2">Pacing amplitude</th><th colspan="2">Projected service life in years (Evera MRI S DR SureScan DDMC3D1)</th></tr><tr><th>Pacing impedance 500 Ω</th><th>Pacing impedance 600 Ω</th></tr><tr><td rowspan="2">DDD, 0%</td><td>2.5 V</td><td>9.7</td><td>9.7</td></tr><tr><td>3.5 V</td><td>9.6</td><td>9.6</td></tr><tr><td rowspan="2">DDD, 15%</td><td>2.5 V</td><td>9.1</td><td>9.2</td></tr><tr><td>3.5 V</td><td>8.7</td><td>8.8</td></tr><tr><td rowspan="2">DDD, 50%</td><td>2.5 V</td><td>8.3</td><td>8.5</td></tr><tr><td>3.5 V</td><td>7.1</td><td>7.4</td></tr><tr><td rowspan="2">DDD, 100%</td><td>2.5 V</td><td>7.3</td><td>7.6</td></tr><tr><td>3.5 V</td><td>5.7</td><td>6.1</td></tr></table>	Pacing mode, percent pacing	Pacing amplitude	Projected service life in years (Evera MRI S DR SureScan DDMC3D1)		Pacing impedance 500 Ω	Pacing impedance 600 Ω	DDD, 0%	2.5 V	9.7	9.7	3.5 V	9.6	9.6	DDD, 15%	2.5 V	9.1	9.2	3.5 V	8.7	8.8	DDD, 50%	2.5 V	8.3	8.5	3.5 V	7.1	7.4	DDD, 100%	2.5 V	7.3	7.6	3.5 V	5.7	6.1
Pacing mode, percent pacing	Pacing amplitude			Projected service life in years (Evera MRI S DR SureScan DDMC3D1)																															
		Pacing impedance 500 Ω	Pacing impedance 600 Ω																																
DDD, 0%	2.5 V	9.7	9.7																																
	3.5 V	9.6	9.6																																
DDD, 15%	2.5 V	9.1	9.2																																
	3.5 V	8.7	8.8																																
DDD, 50%	2.5 V	8.3	8.5																																
	3.5 V	7.1	7.4																																
DDD, 100%	2.5 V	7.3	7.6																																
	3.5 V	5.7	6.1																																
18. Operating conditions	Used in medical and preventive treatment facility. Should be used only by or under the supervision of physicians trained. internal body environment 32 °C-42 °C, humidity 80-100%.																																		
19. Conditions of transportation and storage	Carefully observe these guidelines when handling or storing the device. There are no other requirements for storage and transportation, except for temperature: Store and transport the package between -18 °C and +55 °C.																																		

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	Device handling Checking and opening the package – Before opening the sterile package tray, visually check for any signs of damage that might invalidate the sterility of the package contents. Damaged package – The device packaging consists of an outer tray and an inner tray. Do not use the device or accessories if the outer packaging tray is wet, punctured, opened, or damaged. Return the device to Medtronic because the integrity of the sterile packaging or the device functionality may be compromised. Return the device to Medtronic. This device is not intended to be resterilized. Sterilization These devices are single use only and are not intended to be resterilized. Device temperature – Allow the device to reach room temperature before it is programmed or implanted. Device temperature above or below room temperature may affect initial device function. Dropped device – Do not implant the device if it is dropped on a hard surface from a height of 30 cm or more after it is removed from its packaging. Fluid immersion – Do not immerse the device in fluid or flush the connector ports at the time of implant. Doing so could adversely affect the performance of the device and lead system. “Use by” date – Do not implant the device after the “Use by” date because the battery longevity could be reduced. For single use only – Do not resterilize and reimplant an explanted device. Device storage Avoid magnets – To avoid damaging the device, store the device in a clean area away from magnets, kits containing magnets, and any sources of electromagnetic interference. Temperature limits – Store and transport the package between –18 °C and +55 °C. Electrical reset may occur at temperatures below –18 °C. Device longevity may decrease and performance may be affected at temperatures above +55 °C.
20. Assembly and installation requirements	Not applicable.
21. Warranty	Guaranteed shelf life – 18 months Limited warranty Manufacturer’s Seven (7) Year Limited Warranty and General Warning provided by Medtronic, Inc. (or such other legal entity as may be referred to as the manufacturer on the labeling of this device) (“Medtronic”) Implantable Cardioverter Defibrillator (the “Device”) limited warranty (the “Limited Warranty”) A. What the Limited Warranty provided by Medtronic covers (1) This Limited Warranty is given by Medtronic and provides the assurances set out below to and for the sole benefit of the purchaser of the Device (the “Purchaser”) for a period of seven (7) years commencing on the date of the original implant (the “Implantation Date”) of the Device in a patient. This Limited Warranty gives the Purchaser specific legal rights. The Purchaser may also have other rights, which vary from jurisdiction to jurisdiction and this Limited Warranty does not affect those rights. In this Limited Warranty, “Purchaser” means either (but never both of) A(2) or A(3): (2) The patient in whom the Device has been implanted if he or she purchased the Device and the Replacement (as defined below) or if the Device and the Replacement were purchased on the patient’s behalf by the implanting hospital, other clinical facility or healthcare professional that implanted the Device or the Replacement (“Healthcare Institution”), provided in such a case that the cost of the Device and the Replacement has in each case been charged to the patient’s account by the Healthcare Institution at the Implantation Date and at the date of the implantation of the Replacement, including in circumstances where the patient’s account

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

has been, and/or will be, discharged either wholly or in part by a third party payor (this category of Purchaser is referred to in this Limited Warranty as a "Patient-Purchaser"); or (3) The Healthcare Institution that purchased and implanted the Device at the Implantation Date, provided that:

- a) the cost of the Device has been charged to the account of the relevant Healthcare Institution; and
- b) the Healthcare Institution's claim under this Limited Warranty is solely for the purpose of obtaining a second Medtronic device to be used by it as a Replacement (as defined below) for the Device for the same patient (this category of Purchaser is referred to in this Limited Warranty as an "Institution-Purchaser").

The assurances and benefits to which the Purchaser is entitled under this Limited Warranty are strictly limited to, and do not extend beyond, one or other of the warranty credits set out at A(4) and A(5) below:

(4) Should the Device function in a manner inconsistent with its intended operation and performance due to the quality of materials or workmanship, and should such event occur within the four (4) year period commencing on the Implantation Date, Medtronic will, at its option, either:

- a) Issue a credit to the Purchaser for a Medtronic device to be used as a replacement of the Device (the "Replacement") for the same patient, equal to the lesser of the net price paid by the Purchaser: (i) for the Device (the "Original Purchase Price"), or (ii) for the Replacement (the "Actual Purchase Price"); or,
- b) Without charge, provide to the Purchaser a Replacement for use in the patient in whom the Device was originally implanted that is functionally comparable to the Device.

(5) Should the Device function in a manner inconsistent with its intended operation and performance due to the quality of materials or workmanship, and should such event occur after four (4) years but before the end of seven (7) years from the Implantation Date, Medtronic will provide the Purchaser with a credit against the purchase of the Replacement, in an amount equal to one-half (1/2) of the lesser of the Original Purchase Price or the Actual Purchase Price, reduced on a pro rata basis (at the rate of 1/36th per month) over this three (3) year period.

In no event will any warranty credit issued under this Limited Warranty exceed the Original Purchase Price or the Actual Purchase Price.

(6) THIS LIMITED WARRANTY IS NOT A REPRESENTATION THAT THE BATTERY OF THIS DEVICE WILL LAST THE ENTIRE SEVEN (7) YEAR WARRANTY PERIOD. The Device batteries have a specified capacity that may deplete at different rates depending on individual patient Device settings and individual patient requirements for pacing, cardioversion, defibrillation or other Device functions. Although battery cell depletion will occur with time and is not considered to be caused by the quality of materials or workmanship, the warranty credits set out in Subsections A(4) and A(5) above will still apply to a Device that must be replaced during the applicable warranty period due to battery depletion.

B. Instructions on how to claim under this Limited Warranty

For specific instructions about how to obtain a warranty credit for a replaced Device under this Limited Warranty, Purchasers should contact Medtronic at one of the addresses at the end of this warranty document and applicable to the country where your Device or Replacement was implanted, or contact Medtronic at www.medtronic.com. In general, in order to claim under this Limited Warranty, all of the following conditions must be met:

- (1) The Device must be implanted on or before its "Use By" or "Use Before" date in conjunction with Medtronic leads or leads of equal quality and comparable electrical characteristics.
- (2) The Device must be used in accordance with the labeling and not altered or subjected to misuse, reuse, abuse, accident, or improper handling.
- (3) In jurisdictions where it is a legal requirement or established practice to register a Device with the manufacturer, all Device

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

registration materials must be completed and returned to Medtronic within thirty (30) days of the Implantation Date.

(4) Replaced Devices must be returned to Medtronic at the address listed below within ninety (90) days of explantation. By returning the replaced Device and seeking a warranty credit under this Limited Warranty, the Patient-Purchaser agrees that the replaced Device shall be the property of Medtronic. For Institution-Purchasers, their return of the replaced Device and seeking of a warranty credit under this Limited Warranty is confirmation that they have obtained the prior [written] agreement of the patient that the replaced Device shall be the property of Medtronic.

(5) Replaced Devices must be accompanied by a written statement from the Purchaser indicating that the replaced Device is being returned to Medtronic to determine whether a warranty credit is due under the Limited Warranty.

(6) Where the Purchaser is established or resident in a country within the European Economic Area, the Device should have been first placed on the market in the European Economic Area before a jurisdiction outside the European Economic Area, by Medtronic or by an entity in the same corporate group as Medtronic.

(7) The remedies provided under this Limited Warranty are offered for the sole benefit of the Purchaser when all the above conditions specified in this Limited Warranty are met. The determination that a component of the Device is not as warranted in this Limited Warranty shall be made solely by Medtronic after its tests and inspections and Medtronic's decision shall be final.

THE PURCHASER SHALL ENSURE THAT IT HAS COMPLIED WITH THESE CONDITIONS IN ORDER TO ENSURE THE AVAILABILITY OF THIS LIMITED WARRANTY.

C. This Limited Warranty is limited to its express terms. In particular:

(1) This Limited Warranty is made solely for the benefit of the Purchaser.

(2) The warranty credits provided under this Limited Warranty as described in subsections A(4) or A(5) above are the Purchaser's sole entitlement under this Limited Warranty and no such warranty credit shall extend beyond the period specified in subsections A(4), A(5), and A(6).

(3) To the maximum extent permitted by law, the warranty credits set out at A(4) and A(5) above under this Limited Warranty are exclusive and in lieu of all other warranties and remedies, and Medtronic specifically disclaims all statutory or implied warranties or conditions, including but not limited to, warranties of merchantability, satisfactory quality and fitness for a particular purpose. If statutory or implied warranties cannot lawfully be disclaimed or restricted in a particular jurisdiction, then to the extent permitted by law in such jurisdiction, all such warranties shall be limited in duration to the longer of the legal minimum period or the periods specified in subsections A(4), A(5), and A(6) of this Limited Warranty and to the replacement and credit provisions detailed in subsections A(4) and A(5) of this Limited Warranty.

(4) To the maximum extent permitted by law, Medtronic is not responsible for, and this Limited Warranty is not intended to confer, any responsibility for any direct, special, incidental or consequential damages, losses or expenses resulting from any breach of this Limited Warranty or under any other legal theory whether or not Medtronic was advised or aware of the possibility of such damage, losses or expenses. In some jurisdictions, the foregoing limitation is not legally permitted to apply to death or personal injury claims, fraud or any statutory liability for intentional and grossly negligent or negligent acts and/or omissions, so the above exclusion or limitation may not be applicable to the Purchaser in such jurisdictions.

(5) If any part or term of this Limited Warranty is held to be illegal, unenforceable, or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of the Limited Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Limited Warranty did not contain the particular part or term held to be invalid.

(6) No person has any authority to bind Medtronic to any representation, condition, or warranty with respect to the Device that deviates in any respect from this Limited Warranty.

(7) THIS LIMITED WARRANTY IS NOT APPLICABLE TO LEADS OR ACCESSORIES USED WITH THE DEVICE.

Not to be reproduced without written permission from Medtronic, Inc.

Medtronic, Inc.

- Confidential

	(8) Any warranty issued by the local Medtronic office in the country of purchase supersedes this Limited Warranty. General Warning Medtronic devices, including but not limited to the Device are implanted in the extremely hostile environment of the human body. In addition, each patient requires individual treatment by the physician including but not limited to individual settings of the Device. The hostile body environment and the clinical necessity for patient individual settings for the Device place severe limitations on the design and function of the Device and the lead. These limitations unavoidably reduce the potential performance and longevity of the Device and lead despite the exercise of due care in design, component selection, manufacture, and testing prior to sale. Reference is hereby made to published data on predictable failure rates of devices and the lead and their implantation, which have either been furnished by Medtronic to physicians or are otherwise available to physicians. Important Note Medtronic recommends physicians fully disclose all safety and performance risks associated with the implant of a device or the lead to the patient before the implantation date. Medtronic similarly recommends patients always talk with their physicians about diagnosis and treatment information and ensure that they understand and carefully follow that information. In addition, Medtronic refers physicians to the implant instructions delivered with every device and lead which include further important measures which physicians should take regarding implant techniques and selection of appropriate individual patient settings for a device in order to minimize the risk of reduced safety or performance, including longevity, of a device and lead.																											
22. The list of international regulatory documents/standards, which matches medical product	<table><tr><th>REFERENCED STANDARDS</th><th>TITLE</th><th>HARMONIZATION STATEMENT</th><th>COMPLIANCE</th></tr><tr><td>EN 45502-1:2015</td><td>Active implantable medical devices - Part 1: General requirements for safety, marking and for information to be provided by the manufacturer</td><td>Harmonized under the AIMD and MD Directives of Europe</td><td>Full Compliance</td></tr><tr><td>EN 45502-2-1: 2003</td><td>Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia</td><td>Harmonized under the AIMD and MD Directives of Europe</td><td>Only Section 6 applies</td></tr><tr><td>EN 45502-2-2: 2008</td><td>Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)</td><td>Harmonized under the AIMD and MD Directives of Europe</td><td>Full Compliance</td></tr><tr><td>ISO 14708-1: 2012</td><td>AIMD- Part 1: General Requirements for Safety, Marking and Information to be provided by the Manufacturer</td><td>—</td><td>Full Compliance</td></tr><tr><td>ISO 14708-2: 2012</td><td>Implants for surgery - Active implantable medical devices - Part 2: Cardiac pacemakers</td><td>—</td><td>Only Section 6 applies</td></tr></table>				REFERENCED STANDARDS	TITLE	HARMONIZATION STATEMENT	COMPLIANCE	EN 45502-1:2015	Active implantable medical devices - Part 1: General requirements for safety, marking and for information to be provided by the manufacturer	Harmonized under the AIMD and MD Directives of Europe	Full Compliance	EN 45502-2-1: 2003	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia	Harmonized under the AIMD and MD Directives of Europe	Only Section 6 applies	EN 45502-2-2: 2008	Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance	ISO 14708-1: 2012	AIMD- Part 1: General Requirements for Safety, Marking and Information to be provided by the Manufacturer	—	Full Compliance	ISO 14708-2: 2012	Implants for surgery - Active implantable medical devices - Part 2: Cardiac pacemakers	—	Only Section 6 applies
REFERENCED STANDARDS	TITLE	HARMONIZATION STATEMENT	COMPLIANCE																									
EN 45502-1:2015	Active implantable medical devices - Part 1: General requirements for safety, marking and for information to be provided by the manufacturer	Harmonized under the AIMD and MD Directives of Europe	Full Compliance																									
EN 45502-2-1: 2003	Active implantable medical devices - Part 2-1: Particular requirements for active implantable medical devices intended to treat bradyarrhythmia	Harmonized under the AIMD and MD Directives of Europe	Only Section 6 applies																									
EN 45502-2-2: 2008	Active implantable medical devices - Part 2-2: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (includes implantable defibrillators)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance																									
ISO 14708-1: 2012	AIMD- Part 1: General Requirements for Safety, Marking and Information to be provided by the Manufacturer	—	Full Compliance																									
ISO 14708-2: 2012	Implants for surgery - Active implantable medical devices - Part 2: Cardiac pacemakers	—	Only Section 6 applies																									

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

ISO 14708-6: 2010	Implants for surgery - Active implantable medical devices - Part 6: Particular requirements for active implantable medical devices intended to treat tachyarrhythmia (including implantable defibrillators)	—	Full Compliance
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices (ISO 14971: 2007, corrected version 2007-10-01)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
ISO 11607-1:2009	Packaging for terminally sterilized medical devices - Part 1: Requirements for materials, sterile barrier systems and packaging systems (ISO 11607-1:2006)	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
ISO 11607-2:2006	Packaging for terminally sterilized medical devices Part 2: Validation requirements for forming, sealing and assembly processes - First edition	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN 1041:2008+2013	Information supplied by the manufacturer with medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 11135-1:2014	Sterilization of medical devices – Validation and routine control of Ethylene Oxide sterilization	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN 980:2008	Graphical symbols for use in the labelling of medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-1:2009/AC:2010	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process	Harmonized under the AIMD and MD Directives of Europe	Full compliance for applicable sections: (Meets all requirements for implanted devices in tissue/bone listed in tables 1 and 2)
EN ISO 10993-2:2006	Biological evaluation of medical devices - Part 2: Animal welfare requirements		Full compliance

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

EN ISO 10993-3:2014	Biological evaluation of medical devices - Part 3: Tests for genotoxicity, carcinogenicity and reproductive toxicity	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-4:2017	Biological evaluation of medical devices - Part 4: Selection of tests for interactions with blood	Harmonized under the AIMD and MD Directives of Europe	Not available
EN ISO 10993-5:2009	Biological evaluation of medical devices - Part 5: Tests for in vitro cytotoxicity	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-6:2016	Biological evaluation of medical devices - Part 6: Tests for local effects after implantation	Harmonized under the AIMD and MD Directives of Europe	Full compliance
EN ISO 10993-7:2008/AC:2009	Biological evaluation of medical devices - Part 7: Ethylene Oxide sterilization residuals	Harmonized under the AIMD and MD Directives of Europe	Full Compliance for applicable sections (section 4.3.2 Permanent contact device)
EN ISO 10993-9:2009	Biological evaluation of medical devices - Part 9: Framework for identification and quantification of potential degradation products	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-10:2010	Biological evaluation of medical devices - Part 10: Tests for irritation and skin sensitization		Full Compliance
EN ISO 10993-11:2006	Biological evaluation of medical devices - Part 11: Tests for systemic toxicity	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-12:2012	Biological evaluation of medical devices - Part 12: Sample preparation and reference materials	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-13:2010	Biological evaluation of medical devices - Part 13: Identification and quantification of degradation products from polymeric medical devices	Harmonized under the AIMD and MD Directives of Europe	Full Compliance
EN ISO 10993-15:2000	Biological evaluation of medical devices. Part 15. Identification and quantification of degradation products from metals and alloys	Harmonized under the AIMD and MD Directives of Europe	Full Compliance

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

EN ISO 10993-16:2010	Biological evaluation of medical devices - Part 16: Toxicokinetic study design for degradation products and leachables	Harmonized under the AIMD and MD Directives of Europe	The evaluation was not carried out due to the availability of data on the use of these materials in marketable products and biosafety testing
EN ISO 10993-17:2002	Biological evaluation of medical devices. Part 17. Establishment of allowable limits for leachable substances	Harmonized under the AIMD and MD Directives of Europe	The evaluation was not carried out due to the availability of data on the use of these materials in marketable products and biosafety testing
EN 556-1: 2001/AC:2006	Sterilization of medical devices - Requirements for medical devices to be labeled "Sterile" - Part 1: requirements for terminally sterilized medical devices	—	Full Compliance
ISO 5841-3: 2013	Implants for surgery - Cardiac pacemakers - Part 3: Low-profile connectors (IS-1) for implantable pacemakers	—	Full compliance for applicable sections (Section 4.2 is applicable to leads only and does not apply)
ISO 11318	Cardiac defibrillators – Connector assembly DF-1 for implantable defibrillators -- Dimensions and test requirements	—	Full Compliance
EN 62366-1:2015	Medical devices - Application of usability engineering to medical devices	Harmonized under the MD Directives of Europe	Full compliance only for user interface and other design outputs developed on the Evera MRI project
EN 62304:2016	Medical device software. Software life-cycle processes	Harmonized under the AIMD and MD Directives of Europe	Full compliance only for code developed on the Evera MRI project
EN 60601-1: 2006 AC:2010	General requirements for safety of Medical Electrical devices	Harmonized under the AIMD and MD Directives of Europe (1st Edition)	Clause 14 only (programmer software)

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

	EN 60601-1-6: 2010	Medical electrical equipment -- Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability	Harmonized under the AIMD and MD Directives of Europe	Full compliance only for user interface and other design outputs developed on the Evera MRI project
23. Customer Service	Service Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed below. Medtronic Inc., USA Information for patients: Toll-free +1 800 551 5544 (7 a.m. – 6 p.m., Monday–Friday, central time) Fax +1 763 514 1855 Representative company of Medtronic in Russia Medtronic LLC, Russia, 123317 Moscow, Presnenskaya Naberezhnaya, 10 Tel.: +7 495 580 7377 +7 495 580 7377 FREE Fax: + 7 495 580 7378 E-mail: info.russia@medtronic.ru Internet: www.medtronic.eu, www.medtronic.ru To find the addresses of offices worldwide Medtronic, please go to the following link: http://www.medtronic.com			

Not to be reproduced without written permission from Medtronic, Inc.
Medtronic, Inc.
- Confidential

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www.medtronic.com/manuals

Medtronic, Inc., 710 Medtronic Parkway
Minneapolis, MN 55417 USA

Manufactured in:
Tolochenaz, Switzerland
www.medtronic.com

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


MRP0000AC01 A - L045128AC01 A

2000-00-00
YYMM DD

Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W
Medtronic DDMC3D1		Evera MRI™ S DR SureScan™ [SN] WWW000000W

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