Kimberly J. Horsch
January 9, 2023

"APPROVE"

Medtronic Inc., USA

Pr Regulatory Affairs Specialist

(position)

Ashley Luloff

(name)

Ashley Luloff

(signature)

« 09 » 01 2023

«day» month (numerals)

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2) Foramen needle InterStim (8.8 and 12.6 cm)

The Foramen needle is a unipolar needle that provides intraoperative stimulation of the sacral nerve on the implant. The needle is available in two lengths: 8.8 cm (also referred to as 3.5 inches) and 12.6 cm (also referred to as 5 inches). The figure below provides information about the design of the needle


Figure 7. Foramen needle with pre-loaded stylet.

- ① Needle heel
- ② Needle heel (side view)
- ③ Non-insulated section
- ④ Black dot
- ⑤ Needle hub
- ⑥ Needle stylet hub

Figure 3 - Foramen needle InterStim

3) Test stimulation cable, model 3575

The test stimulation cable should be used during the acute test stimulation phase to connect an external neurostimulator to the foramen needle.

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


Figure 4 - Test stimulation cable

4) Ground pad, model 041826

Ground pads are silicone rubber electrodes that self-adhere to the patient's skin using a conductive hydrogel. After application, the pads form a positive lead when connected to an external neurostimulator.

Non-sterile, self-adhesive ground pad, used during lead implantation and/or test stimulation procedures. Connects to test stimulation cable and patient cable.

The non-sterile ground pad is for one patient only, but it can be reapplied to the same patient multiple times. The clinician or patient may, as appropriate, reapply a new ground pad during a temporary sacral nerve stimulation procedure in the clinical conditions or during an ongoing test at home.

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
Figure 5 – Ground pad

5) Test stimulation lead, model 3059

The InterStim Test Stimulation Lead is used for temporary sacral neuromodulation.

The Test Stimulation Lead is a single-use temporary device for one patient, intended to be implanted for a maximum of seven (7) days.

Figure 6 – InterStim Test Stimulation Lead

6) Quadripolar lead 28/33/41 cm for electrical stimulation, model 3889-28/33/41

The quadripolar lead 28/33/41 are quadripolar in-line leads designed to be implanted adjacent to the sacral nerve. The lead contains four cylindrical electrodes equal in length and spaced equidistantly.

The leads have tines and marker bands. The tines anchor the lead, and the marker bands indicate lead depth and tine deployment during percutaneous implantation with a Medtronic lead introducer. Marker bands C and D can be seen visually, while marker bands A and B can be seen only under fluoroscopy.

The leads have colored marker tips at the proximal end to ensure the lead is fully seated in the Neurostimulator connector block.

The InterStim system, which is used with a quadripolar electrical stimulation lead, is only suitable for head MRI scans.

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

Figure 7 – Quadripolar lead for electrical stimulation


7) InterStim SureScan MRI lead kit – electrode length 28/33/41 cm

The InterStim SureScan MRI leads are quadripolar inline leads designed to be implanted adjacent to the sacral nerve. The electrodes are equal in length and spaced equidistantly. The leads have tines and marker bands. The tines anchor the lead, and the marker bands indicate lead depth and tine deployment during percutaneous implantation with a Medtronic lead introducer. Marker bands C and D can be seen visually, while marker bands A and B are intended for visualization under fluoroscopy.

The leads are blue at their proximal ends, which serve as visual aids to help determine that they are fully seated when connecting to either the percutaneous extension or the neurostimulator.

Unlike standard leads (see item above), shielded MRI-compatible leads (InterStim SureScan MRI lead, model 978B1) reduce the risk of thermal tissue damage by redistributing RF energy along the entire length of the lead.

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- | | |
|--|-----------------------------|
| ① Bent-tip stylet directional
indicator | ⑥ Electrode numbers |
| ② 4.32 mm spacing (proximal
end) | ⑦ Connector numbers |
| ③ Removable mini-hook hook
guide | ⑧ Marker band B |
| ④ Marker band D | ⑨ Marker band A |
| ⑤ Marker band C | ⑩ 3 mm spacing (distal end) |

Figure 8 – InterStim SureScan MRI lead


Figure 9 – Stylet

The lead kits contain the following:

- Quadripolar lead (with inserted bent-tip stylet and removable mini-hook guide)
- Straight stylet
- Introducer with dilator
- Tunneling tool, tunneling tip, and tube
- Ground pad (nonsterile)
- Torque wrench
- Test stimulation cable with mini-hook
- Directional guide
- Two 3.5-inch foramen needles with stylets (5-inch foramen needles packaged and sold separately)


The contents of the inner package are sterile. The lead and all accessories are intended for single use only.

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The percutaneous extensions are used for the advanced evaluation of the sacral neuromodulation system. The percutaneous extension is for use with the InterStim SureScan MRI lead.

Each extension has a connector at its distal end that connects to a Medtronic neurostimulation lead, and a connector at its proximal end that connects to a external neurostimulator (ENS).


Figure 10 – InterStim Percutaneous Extension

9) Lead Introducer, model 042294

The Medtronic Model 042294 Lead Introducer and Model 3550-18 Lead Introducer Kit are used to introduce compatible Medtronic Neurological leads used for Sacral Nerve Stimulation (SNS).

The Model 042294 Lead Introducer contains six lead introducers. Each lead introducer contains one directional guide, one dilator, and one introducer sheath.

The Model 3550-18 Lead Introducer Kit contains:

- Lead introducer – 1 pc.;
- Foramen needle, 9 cm – 2 pc.;
- Ground pad – maximum 2 pc.;
- Long test stimulation cable – 1 pc.;
- Patient cable – 1 pc.

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
Figure 11 – Introducer sheath (top), directional guide (center), dilator (bottom)

10) Twist lock cable 25, 64 cm Model 357625/357664

The temporary percutaneous extension connects the permanent lead to the external neurostimulator using a twist lock cable.


Figure 12 – Twist lock cable

11) Percutaneous Extension and Tunneling Tools Accessory Kit model 3550-05

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


Figure 13 – Percutaneous Extension and Tunneling Tools Accessory Kit

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The Percutaneous Extension and Tunneling Tools Accessory are designed for use with compatible Medtronic lead kits and are for single use only. The kit includes a linear extension, a tunneling guidewire, a guidewire tube and one tip.

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The Percutaneous Extension is used for extended evaluation of sacral neuromodulation therapy. The extension is used with a quadropolar lead 28/33/41 cm for electrical stimulation. Each extension at the distal end has a connector that connects to a Medtronic neurostimulation lead, and the connector at the proximal end connects to a twist lock cable, which in turn connects to an external neurostimulator (ENS).

12) iCon patient programmer

The iCon patient programmer can be used to control the InterStim II neurostimulator. Patient programmers are portable, battery-powered devices that patients use to control the stimulation of an implanted neurostimulator. The patient programmer allows patient to:

- Turn your neurostimulator on and off.
- Adjust amplitude within the limits set by the clinician.
- Change programs within the programs set by the clinician.
- Warn when the battery of the neurostimulator is low.
- Provide patient profile and therapy usage information (only available with a Model 3037 Patient Console that is linked to a neurostimulator)

Figure 16 – iCon patient programmer

The programmer interface does not have access to control settings that can only be set by a physician (eg pulse width, pulse frequency, lead configuration, etc.).

The programmer allows the patient to:

- turn on the programmer by pressing the on/off button. stimulation or increase/decrease buttons;
- immediately turn off the stimulation by touching the on/off button. pacing and selecting the pacing off option that appears on the touch screen.
- increase and decrease stimulation levels;

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- track the battery levels of the Neurostimulator and programmer.
- switch from one available stimulation mode set by the doctor to another (up to four), in accordance with his instructions.
- "lock" the programmer - this function actually cuts off the power supply of the programmer and inactivates the touch screen
- View troubleshooting messages that inform the patient of problems that may interrupt stimulation, including:
 - stimulation off
 - the batteries of the Neurostimulator are low
 - low charge of the Neurostimulator
 - stuck button
 - Batteries of the programmer are low
 - low batteries of the programmer
 - check connection: notification to the patient that the cable has become disconnected
 - system malfunctions

After synchronizing, the Therapy screen appears

Therapy screen icons:

Icons	Description
	Neurostimulator is on
	Neurostimulator is off
	Synchronizes the neurostimulator and patient programmer
	Decreases amplitude.
	Increases amplitude.

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HH90 Handset

The patient Mobile Platform Handset is a non-medical, commercially available, off-the-shelf Samsung handset with Android operating system. The features and software available on the patient Mobile Platform Handset will be controlled by mobile device management (MDM) software. The Samsung handsets use Bluetooth Low Energy (BLE) 4.0 for communication with the Model TM90 Communicator. This is a non-medical device.

TM90 Communicator

The communicator acts as a wireless communication "bridge" between the application installed on the handset and the implantable medical device. The implanted device communicates with the application only through the communicator.

The exchange of information between the application and the implanted device via the communicator can only take place after the so-called pairing of the application with the communicator has been performed. Pairing is done before using the app for the first time.

The Medtronic Model TM90 Communicator is intended for patient or clinician use, to connect the app installed on the Model HH90 handset to the implanted medical device.

The app on the handset interacts with the implanted medical device through the communicator.

The communicator contains a lithium-ion battery that will need to be regularly charged. The communicator can be used in a home or health care setting.


Figure 17 - Communicator

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The communicator includes a power button and USB charging port. A battery indicator light shows when the battery is low (yellow), charging (orange), or fully charged (green).

The BLUETOOTH® technology indicator light (blue) shows if the communicator is connected to the handset.

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Продолжительность циклов включения/выключения ^c	0,1 с	от 0,1 с до 24 ч	0,1 с	Output On/Off	Off	On or Off	Off
Функция SoftStart/Stop ^d	Выкл.	1, 2, 4, 8, 15, 30 с или Выкл.	Выкл.	Polarity			
Выход Вкл./Выкл.	Выкл.	Вкл. или Выкл.	Выкл.	Electrode Number			
Полярность Номер электрода				0		Off, + or –	
				1	Off	Off, + or –	Off
				2	Off	Off, + or –	Off
				3	+	Off, + or –	+
				Case ^e	Off	Off or +	Off
0		Выкл., + или –		^a All values are approximate. ^b Power-on-reset (POR) turns OFF stimulation by resetting the amplitude to 0.0 V. POR can occur when there is a temporary fluctuation in battery voltage (eg, due to electromagnetic interference during electrocautery or defibrillation) or the battery is depleted. When POR occurs, the serial number is reset to a nominal value and must be entered with the clinician programmer. ^c Cycle On and Cycle Off time must be greater than or equal to the SoftStart/Stop time. ^d SoftStart/Stop must be less than or equal to the Cycle On and Cycle Off time. ^e If case electrode is positive, electrodes 0-3 may be programmed to Off or – only.			
1	Выкл.	Выкл., + или –	Выкл.				
2	Выкл.	Выкл., + или –	Выкл.				
3	+	Выкл., + или –	+				
Корпус ^e	Выкл.	Выкл. или +	Выкл.				

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Table 2. Physical characteristics of the InterStim II Neurostimulator

Description	Value
Height	44 mm (1.7 in)
Length	51 mm (2 in)
Thickness	7.7 mm (0.3 in)
Weight	22 g (0.77 oz)
Volume	12.5 cm ³ (0.76 in ³)
Surface area	47.3 cm ² (7.33 in ²)
Power source	1.3 Amp Hours, 3.2 V Lithium silver vanadium oxide hybrid
Expected lifetime	4.8 to 6.3 years

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
Table 4 - Technical characteristics of the foramen needle

Parameter	Value
М (Foramen needle 8.8 cm InterStim)	88 mm
М (Foramen needle 12.6 cm InterStim)	126 mm
K_{max} Foramen needle 8.8 cm InterStim)	116,6 mm
K_{max} (Foramen needle 12.6 cm InterStim)	167,4 mm
Needle sharpening angle along the main plane	$12^\circ \pm 2^\circ$
Requirements for the outer surface of the needle	The outer surface of the needle must be smooth, light, uniform, without cracks, shells, nicks, scratches, chipped places, delaminations, burrs and other visible defects. The stylet (mandrin) should freely enter the needle channel.
Straightness deviation	at least 0,2 mm at length 50 mm
Roughness parameter R_a of outer surface	at least 0,16 μm
Leakage current through the insulation of the conductor needle	at least 1 μA

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4 Test stimulation cable, model 3575


Table 5 - Specifications of the test stimulation cable

Parameter	Maximum value, mm
Cable length (from ENS connector to hook)	$2286 \pm 76,2$
Length of the contact part (hook side)	381
Connector length	$45 \pm 0,8$
Length of the contact part (from the side of the connector for ground electrode)	$444,5^{+50,8}_{-25,4}$
Cable length (from ENS connector to earth electrode connector)	$698,5 \pm 25,4$

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5 Ground pad, model 041826


Table 6 - Technical characteristics of the ground pad

Parameter	Value, mm
Length	85,5±3
Width	57,0±3
Pad length	38,0±1,5
Pad width	44,5±1,5

Ground pad resistance, not less than - 150 Ohm

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6 Test stimulation lead, model 3059


Table 7 - Technical characteristics of the electrode for test stimulation

Parameter	Maximum value, mm
Lead length	410
The length of the lead contact part	10
Lead diameter	0,66

Test stimulation lead resistance, not less than - 150 Ohm

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7 . Quadripolar lead 28/33/41 cm for electrical stimulation, model 3889-28/33/41

The total length of the lead, depending on the model:

280±10% mm,

330±10% mm,

410±10% mm

Weight of the quadripolar electrode for electrical stimulation:

Not more than 2 g

The quadripolar lead for electrical stimulation withstands the following loads without destruction of any of the conductors or the formation of cracks in the electrical insulation or in the body of the conductor:

a) tensile load, not less than 44.0 N.

For verification, the electrode cable is loaded with the specified load 25 times in the direction along the axis. The load is applied without jerks for 1 s each time.




b) bending load, not less than 25.0 N.

To check, a mass equal to 10D2 grams is suspended from the free end of the cable, where D is the diameter of a round cable, for a flat cable - the minimum cable size, mm.

The electrode cable is bent at an angle of 90°±10° in the plane of least resistance. Bending is done 25 times in each direction.

c) the maximum fixing force of the electrode connection with the Product, according to ISO 14708-1, is not less than 7 N.

To check, the assembled pair of electrode with the Product is immersed in a saline solution with a concentration of about 9 g/l at room temperature for at least 8 hours. 1 s. In this case, the pair should not be disconnected.

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Table 9 - Technical characteristics of the Tunneling tool, tunneling tip, and tube

Parameter	Value
Tunneling tool length:	206,8 mm $\pm$ 0,8 mm
Tube length:	198,6 m $\pm$ 1,0 mm
Tube diameter:	3,43 mm $\pm$ 0,15 mm
Tunneling tip length:	25,4 mm $\pm$ 1,5 mm
Tunneling tip diameter:	4,0 mm $\pm$ 0,1 mm
Tunneling tip angle:	15° $\pm$ 3°

See 10.10 for technical specifications of the introducer
See 10.5 for technical specifications of the Ground pad
See 10.1 for technical specifications of the torque wrench
See 10.4 for technical specifications of the Test stimulation cable
See 10.3 for technical specifications of the Foramen needle

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9 InterStim Percutaneous Extension


Table 10 – Technical characteristics of the InterStim Percutaneous Extension

Parameter	Value
Length	1000 mm ± 50,8 mm
Extension body diameter	1,32 mm ± 0,076 mm
Distal end (lead side)	
Connector	Quadripolar
Number of electrical conductors	4
Distal connector length	29,2mm ± 0,5mm
Distal connector diameter	6,0 mm ± 0,1 mm
Conductor wire diameter	0.15 mm (0.006 in)
Lead entrance diameter	2.5 mm (0.065 in)
Connector (at the proximal end on the side of the neurostimulator)	Custom external neurostimulator interface, 4 active conductors

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10 Lead Introducer, model 042294 and Lead Introducer kit, model 3550-18


Table 11 – Technical characteristics of the Lead Introducer

Parameter	Value
Introducer Sheath length (not including adapter and handle)	110 mm ±1,3 mm
Introducer Sheath diameter (inner)	1,67 mm ±0,03 mm
Introducer Sheath diameter (outer)	2,49 mm
Dilator length (not including adapter and handle)	124,8 mm ± 2,5 mm
Dilator diameter (inner)	0,69 mm
Dilator diameter (outer)	1,85 mm
Guide wire length	304,8 mm ± 2,5 mm


Table 12 – Technical characteristics of the Patient cable

Parameter	Max Value, mm
Length of cable	1355
Contact part length	58
Connector length	27
Width of the base of the contact part	13,9
Connector Width	9,5

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
Table 13 – Technical characteristics of the Foramen needle, 9 cm

Параметр	Max value, cm
M, model 9 cm	9,0
K _{max} model 9 cm	12,6

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11 Twist lock cable 25/64 cm Model 357625/357664


Table 14 – Technical characteristics of the Twist lock cable

Parameter	Minimum value, mm
Length of the size 25 cm	250
Length of the size 64 cm	640

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12 Percutaneous Extension and Tunneling Tools Accessory Kit model 3550-05


Table 15 – Technical characteristics of the linear extension

Parameter	Value
Extension length	330,2 mm
Extension connector length	21 mm
Connector outer diameter	6,2 mm
Connector inner diameter	1,5 mm
Number of electrical conductors	4

Table 16 – Technical characteristics of the tunneling guidewire, tip, guidewire tube

Parameter	Value
Tunneling tool length:	206,8 mm ± 0,8 mm
Tube length:	198,6 mm ± 1,0 mm
Tube diameter:	3,43 mm ± 0,15 mm
Tunneling tip length:	25,4 mm ± 1,5 mm
Tunneling tip diameter:	4,0 mm ± 0,1 mm
Tunneling tip angle:	15° ± 3°

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13 iCon patient programmer

The programmer has the following overall dimensions:

Table 17 - Overall dimensions of the programmer

Parameter	Value, mm, no more than
Length	94,0
Width	56,0
Height	28,0

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
Table 18 - Overall dimensions of the carrying case

Parameter	Minimum value, mm
Length	100
Width	60
Height	30

- 1) The weight of the programmer is not more than 120 g.
- 2) The programmer is powered from an independent source, voltage 3.0 V, capacity not less than 2.0 Ah - 2 AAA batteries.
- 3) The programmer provides operation after synchronization with the Neurostimulator via a radio channel with a frequency band of 160-190 kHz.
- 4) The programmer is a programmable device, programmed using the NVision Physician Programmer with accessories - FSZ 2010/07875, manufactured by Medtronic, Inc., USA
- 5) The programmer, together with the software (SW), ensures operation in the patient mode with the following parameters:
 - activation and deactivation of stimulation;
 - blocking of the programmer;
 - regulation of the amplitude within the limits set by the doctor;
 - indication of the battery level of the Product and the programmer;
 - switching of stimulation modes (up to 4 modes) set by a doctor;
 - view notifications about system malfunctions.
- 6) The parameters of the programmer, except for those specified for the patient mode, are not available for changing by the patient. The boundary conditions of the fixed inputs and outputs are checked by the software. The impossibility of making unforeseen changes confirms that they are properly controlled.

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15 Handset with communicator, model TM90**15.1 Communicator Model TM90**

Table 19 – Technical characteristics of the communicator

Parameter	Value
Power source	Rechargeable battery, USB 2.0 compliant cable with external power supply
Battery life - Clinician use	2 hours per day with power on
Battery life - Patient use	30 minutes per day with power on
Mode of operation	Continuous
Classification	Internally powered
Length	60,25 mm (2,4 in)
Width	9,35 mm (0,37 in)
Height	121,25 mm (4,8 in)
Weight, including batteries	57 g (2 oz)
Service life	2 years minimum ^a
Operating temperature	5°C to 40°C ^{bc} (41°F to 104°F)
Storage temperature	-25°C to 70°C ^d (-13°F to 158°F)

^a Communicator can continue to be used after the minimum service life has elapsed.

^b A relative humidity range of 15% to 90%, noncondensing, but not requiring a water vapour partial pressure greater than 50 hPa; and an atmospheric pressure range of 700 hPa to 1060 hPa.

^c At the upper range of operating temperature, the communicator may reach 43°C (109°F).

^d 5°C to 35°C at a relative humidity up to 90%, noncondensing; and >35°C to 70°C at a water vapour pressure up to 50 hPa.

Near-field magnetic induction communication

The Medtronic Model TM90 communicator uses near-field magnetic induction communication to interact with compatible medical devices. Successful programming of the compatible medical device requires close proximity between the Model TM90 communicator and the medical device.

Table 20 – Wireless operating characteristics: Near-Field Magnetic Induction Communication for Model TM90

Radio type	Near-field magnetic induction communication
------------	---

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Operating frequency band	175 kHz
Effective radiated power output	H-field strength of <30 dBμA/m at 10 m
Receiving bandwidth	Up to 30 kHz (occupied bandwidth of the transmitter)
Modulation (modulation is model-specific)	On-off keying (OOK) burst
Effective range	Keep the TM90 communicator within 4 cm (about 1.5 inches) of the compatible Medtronic neuromodulation medical device to ensure successful programming.

BLUETOOTH® wireless technology

The Medtronic Model TM90 communicator uses BLUETOOTH wireless technology to exchange data and information with the commercial off-the-shelf handset. BLUETOOTH wireless technology is a wideband transmission system, which uses wideband modulation techniques. The use of the protocol by Medtronic includes forward error correction (FEC) on data packets, frequency hopping spread spectrum (FHSS) that aids in spectral co-existence, and a reliable data stream profile. Interrupted or partial messages are invalid and must be resent. If communication errors are detected, the system will automatically retransmit the interrupted message to guarantee delivery of the message. If the wireless connection does not automatically re-establish itself, the user should reposition the Model TM90 communicator to re-establish communication.

Table 21 – Wireless operating characteristics: BLUETOOTH wireless technology for Model TM90

Radio type	BLUETOOTH 4.0
Operating frequency band	2.4 to 2.4830 GHz ISM
Effective radiated power output	<10 mW e.i.r.p.
Receiving bandwidth	2 MHz per aggregated channels (bandwidth period 0.5) (occupied bandwidth of the transmitter)
Modulation (modulation is model-specific)	Gaussian frequency shift keying (GFSK)

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Effective range

For wireless communication, the handset must remain within 2 meters (6 feet) of the communicator.

Table 22 – Technical characteristics of the charger

Parameter	Value
Rated Input	100-240V - 50-60Hz, 0.15A
Rated Output	5.0V 1.0A ^a
Isolation Mark	Class II equipment
Certification	Agency certified
Charging Cord	Micro USB 2.0 standard compliant

^a Minimum

The USB charging cable and AC adaptor are included in the handset package. Use this cable or another USB 2.0 compliant cable to charge the device.

10.15.2 Handset

Table 23 – Technical characteristics of the handset

Parameter	Value
Manufacturer	Samsung Electronics
System	
Operating system	Android 7.0 or higher
Charging connector type	Micro USB
Connecting external devices via USB	USB 2.0
Overall dimensions (Length x Width x Height), mm	147,3 x 69,6 x 8,4
Weight, g	141
IPX	IP 68
Case (Length x Width x Height (including strap)), mm	152,4 x 119,4 x 45,72 (6.0 x 4.7 x 1.8 in)

Table 24 – Technical characteristics of the USB charger

Параметр	Значение
Rated Input	100-240V - 50-60Hz, 0.15A
Rated Output	5.0V 1.0A ^a
Isolation Mark	Class II equipment
USB-cable	Length - 1000 mm ± 50 mm

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
Ashley Lutloff
Pr. Regulatory Affairs Specialist
Medtronic Inc.

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MA12231A061 Rev A

Printing instructions: *Refer to doc#
A88898 for Printing Instructions*

2021-10-15

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.4: 03-30-2021

Medtronic and the Medtronic logo are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

3058 2021-10-15 English 3

MA12231A061 Rev A

2021-10-15

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.4: 03-30-2021

4 English 3058 2021-10-15

2021-10-15

MA12231A061 Rev A

Table of contents

Introduction 7

Description 7

Intended purpose 7

Package contents 7

Device specifications 8

X-Ray identification 12

Implanted components and MRI scans 12

Implant criteria for full-body MRI scan eligibility 12

When changing components 13

When explanting components 13

Instructions for use: Model 3058 14

Verifying neurostimulator operation 14

Creating a pocket for the Model 3058 Neurostimulator 14

Connecting the lead to the Model 3058 Neurostimulator 15

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The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, patient selection, individualization of treatment, and component disposal. For customers in Japan, the appropriate package insert provides information about safety, contraindications, warnings, precautions, and adverse events.

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USA The clinical summary provides information about the clinical study results for the neurostimulation system.

Introduction

This manual provides instruction on using the Medtronic InterStim II Model 3058 Neurostimulator.

For the MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan, refer to the MRI Guidelines for the InterStim system manual.

Description

The Medtronic InterStim II Model 3058 Neurostimulator is part of a neurostimulation system for sacral neuromodulation therapy.

Intended purpose

The implantable neurostimulator generates electrical pulses and delivers stimulation through one lead as part of a neurostimulation system for sacral neuromodulation therapy.

Package contents

- Neurostimulator
- Torque wrench
- Product literature, including patient identification card

Filename Date Time
UC200xxxxxx EN
4.6 x 6 inches (117 mm x 152 mm)

Medtronic Confidential
ImplantManual.xsl - IPGTemplate.fm
Template version 6.4: 03-30-2021

Device specifications

The implantable neurostimulator (INS) is a programmable device that delivers stimulation through a lead. Programmable parameters include amplitude, pulse width, rate, and cycling. Refer to Table 1 for shipping, operating, and power-on-reset values. Refer to Table 2 for physical characteristics and materials and substances.

8 English 3058 2021-10-15

2021-10-15

MA12231A061 Rev A

Table 1. Shipping, operating, and power-on-reset (POR) values for the Medtronic InterStim II Model 3058 Neurostimulator³

Programmable Parameters	Shipping	Operating	POR ^b
Amplitude			
Normal Resolution	0.0 V	100 mV steps	0.0 V
Upper Limit	0.0 V	8.5 V maximum	0.0 V
Lower Limit	0.0 V	0.0 V minimum	0.0 V
Fine Resolution		50 mV steps	
Upper Limit		6.35 V maximum	
Lower Limit		0.0 V minimum	
Rate	14 Hz	49 values (from 2.1 to 130 Hz)	31 Hz
Pulse width	210 μ s	Increments of 30 μ s steps 450 μ s maximum 60 μ s minimum	210 μ s
Operating Mode	Continuous	Continuous or Cycling	Continuous
Cycle On/Cycle Off time ^c	0.1 sec	0.1 sec to 24 hr	0.1 sec
SoftStart/Stop ^d	Off	1, 2, 4, 8, 15, 30 sec, or Off	Off
Output On/Off	Off	On or Off	Off
Polarity			
Electrode Number			
0		Off, + or –	
1	Off	Off, + or –	Off
2	Off	Off, + or –	Off

Table 1. Shipping, operating, and power-on-reset (POR) values for the Medtronic InterStim II Model 3058 Neurostimulator^a (continued)

Programmable Parameters	Shipping	Operating	POR ^b
3	+	Off, + or –	+
Case ^e	Off	Off or +	Off

^a All values are approximate.

^b Power-on-reset (POR) turns OFF stimulation by resetting the amplitude to 0.0 V. POR can occur when there is a temporary fluctuation in battery voltage (eg, due to electromagnetic interference during electrocautery or defibrillation) or the battery is depleted. When POR occurs, the serial number is reset to a nominal value and must be entered with the clinician programmer.

^c Cycle On and Cycle Off time must be greater than or equal to the SoftStart/Stop time.

^d SoftStart/Stop must be less than or equal to the Cycle On and Cycle Off time.

^e If case electrode is positive, electrodes 0-3 may be programmed to Off or – only.

Table 2. Physical characteristics of the InterStim II Model 3058 Neurostimulator.^a

Description	Value
Height	44 mm (1.7 in)
Length	51 mm (2 in)
Thickness	7.7 mm (0.3 in)
Weight	22 g (0.77 oz)
Volume	12.5 cm ³ (0.76 in ³)
Surface area	47.3 cm ² (7.33 in ²)
Materials and substances to which the patient can be exposed ^{b c d e f}	Titanium, polyurethane, silicone, silicone medical adhesive
Power source ^g	1.3 Amp Hours, 3.2 V Lithium silver vanadium oxide hybrid
Expected lifetime	4.8 to 6.3 years
Serial number ^h	
Radiopaque identification (ID)	NJY
Transmitter and receiver	
Carrier frequency	175 KHz
Output power	≤30 dBμA/m
RF Wireless Technology	Near Field Magnetic Inductive Communication

^a All measurements are approximate.

^b Discuss any allergies or other intolerances related to the materials and substances with the patient before procedure.

^c Tested for category 1A or 1B carcinogenic, mutagenic, or toxic for reproduction (CMR) substances, or endocrine disrupting chemicals (EDC). No known CMRs or EDC were found in the tested materials or substances.

^d The materials are ordered from the most to least amount.

^e Does not contain natural rubber latex.

^f Does not contain nickel.

^g The power source is hermetically sealed within the case.

^h The serial number is the radiopaque ID followed by a number. The clinician programmer displays the entire number beginning with the radiopaque ID.

X-Ray identification

Radiopaque identification shows the manufacturer and neurostimulator model number (Figure 1) using standard x-ray procedures. The code appears as black characters on white background. The Medtronic symbol identifies Medtronic as the manufacturer. For the Medtronic InterStim II Model 3058 Neurostimulator, the designated characters are NJY.


Figure 1. The InterStim II Model 3058 Neurostimulator radiopaque code block.

Implanted components and MRI scans

Implant criteria for full-body MRI scan eligibility

To allow a patient full-body MRI scan eligibility under specific conditions, implant a Medtronic neurostimulation system with SureScan™ MRI technology as follows:

- Use only SureScan MRI neurostimulation system components (for example, leads and neurostimulators).
- Explant any previously abandoned sacral neuromodulation therapy leads or extensions that may be in the patient (that is, leads or extensions, or portions of, that are not connected to a neurostimulator).

Note: Confirm MRI compatibility of any other implanted medical devices. Other implanted medical devices may limit or restrict MRI scans.

- Enter all component model number and implant location information using the clinician programmer (also referred to as the handset).

If the above implant criteria are not met, the patient will not have a neurostimulation system with full-body MRI scan eligibility. MRI scan eligibility will be restricted. For the

MRI conditions and MRI-specific warnings and precautions for conducting an MRI scan, refer to the MRI Guidelines for the InterStim system.

When changing components

When adding, changing, or removing neurostimulators, leads, extensions, and accessories, always enter up-to-date information into the programmer.

When explanting components

Warning: If permanently explanting a neurostimulator, be sure to also explant all leads, extensions, and accessories. Abandoned components may prevent the patient from being allowed MRI scans in the future due to concerns of lead electrode heating that can result in tissue damage.

Instructions for use: Model 3058

Cautions:

- When using sharp instruments near the neurostimulator, be careful to avoid nicking or damaging the case or the connector block. Damaging the neurostimulator may require surgical replacement.
- Do not use saline or other ionic fluids at connections, which could result in a short circuit.

Verifying neurostimulator operation

Before opening the sterile neurostimulator package, use the clinician programmer to interrogate the neurostimulator and verify neurostimulator battery status and current settings.

 **Warning:** Avoid implanting a damaged device. Severe burns could result from exposure to the battery chemicals. Damaged devices may fail to operate properly over the service life of the neurostimulator.

 **Warning:** Implant the device as far away as possible from another active implanted device (for example, pacemaker, defibrillator) to minimize possible interaction between the devices. For information on interactions between multiple devices, consult the manufacturer's labeling for risks associated with the other devices.

 **Caution:** Do not implant a neurostimulator if it was dropped onto a hard surface from a height of 30 cm (12 in) or more, because the neurostimulator may be damaged and fail to operate properly over the service life of the neurostimulator.

 **Caution:** Select a neurostimulator implant location that meets the following criteria:

- Away from bony structures (for example, 3 – 4 cm [1.2 – 1.6 in]) to minimize discomfort at the neurostimulator site
- Away from areas of restriction or pressure to minimize the potential for skin erosion, patient discomfort, or damage to components
- In an area accessible to the patient for proper operation of a patient control device

Creating a pocket for the Model 3058 Neurostimulator

1. Create a subcutaneous pocket for the neurostimulator by blunt dissection to the posterior surface of the muscle. The neurostimulator is placed in the upper buttock area.

Cautions:

- Ensure that the neurostimulator is placed no deeper than 2.5 cm (1 in) below the skin and is parallel to the skin. If the neurostimulator is too deep or is not parallel to the skin, telemetry may be unsuccessful.
- If the patient has any other neurostimulator, the neurostimulators must be separated by a minimum of 20 cm (8 in).

 **Caution:** The neurostimulator is provided sterile. Do not soak the neurostimulator in antibiotic solution. Soaking in antibiotic solution can possibly affect lead connections.

Note: To avoid infection, it is recommended that the neurostimulator implant site be irrigated with antibiotic solution, and that IV antibiotics be administered perioperatively. Do not allow the neurostimulator to come into contact with any nonsterile surface. Do not place on skin. If an infection occurs, it may require surgical removal of the implanted system.

2. Place the neurostimulator in the pocket to assure proper fit and then remove it. Keep the neurostimulator sterile and clean.
3. Tunnel from the lead incision site to the neurostimulator pocket. Refer to the product literature packaged with the lead for detailed tunneling and lead implant instructions.
4. Connect the lead to the Model 3058 Neurostimulator according to the steps in "Connecting the lead to the Model 3058 Neurostimulator" in this manual.

Connecting the lead to the Model 3058 Neurostimulator

 **Caution:** Before connecting components, wipe off any body fluids and dry all connections. Fluids in the connection may result in stimulation at the connection site or intermittent stimulation.

1. Wipe the proximal lead electrode connections with sterile gauze. If necessary, use sterile (United States Pharmacopeia [USP]) water or a nonionic antibiotic solution, then wipe dry.
2. Make sure the connector block receptacle is dry and clean.
3. Confirm that the lead has four electrodes matching the encapsulated diagram on the neurostimulator (Figure 2).
4. Insert the lead into the neurostimulator connector block until fully seated (Figure 2).

 **Caution:** Do not pull the lead body taut when implanted. The lead is available in different lengths. Select a length that allows connection without tension.

 **Caution:** Do not insert the lead into the connector block if the setscrew is not sufficiently retracted. If the setscrew is not retracted, the setscrew may damage the lead, or the lead will not seat fully into the connector block.

Notes: If the lead does not seat fully into the connector block, take the following actions:

- Insert the torque wrench into the self-sealing grommet and rotate the setscrew 1/2 turn counterclockwise (Figure 3).
- Then insert the lead fully into the connector block.
- If unsuccessful, repeat above a maximum of three times. Do not turn the setscrew more than two total turns to prevent the removal of the setscrew from the connector block.


Figure 2. Insert lead fully into Model 3058 Neurostimulator connector block.

① Encapsulated diagram


Figure 3. To retract out the setscrew from the Model 3058 Neurostimulator, use the torque wrench and turn setscrew counterclockwise.

Note: The torque wrench must be oriented to the same angle as the setscrew (Figure 3).

5. Fully insert the torque wrench into the self-sealing grommet of the connector block and tighten the setscrew by turning clockwise until the torque wrench clicks (Figure 4).

Cautions:

- Be sure the torque wrench is fully inserted into the self-sealing grommet. If the torque wrench is not fully inserted, the setscrew may be damaged, resulting in intermittent or loss of stimulation.
- Before tightening the setscrew, ensure that the lead is fully inserted into the connector block to prevent damaging the lead.
- Verify that the self-sealing grommet is closed after the torque wrench is withdrawn. If fluid leaks through a grommet seal that is not fully closed, the patient may experience shocking, burning, or irritation at the neurostimulator implant location, or intermittent stimulation or loss of stimulation may occur.

Note: Discard the torque wrench after surgery. The torque wrench is single-use-only. Its operation cannot be assured if it is used for multiple surgeries.


Figure 4. Tighten the setscrew in the self-sealing grommet by turning clockwise until the torque wrench clicks.

Note: The torque wrench must be oriented to the same angle as the setscrew (Figure 4).

Implanting the Model 3058 Neurostimulator

1. Place the neurostimulator into the subcutaneous pocket with the etched identification side placed outward, away from muscle tissue, and ensure that the lead is not bent sharply and has sufficient strain relief.

 Cautions:

- Ensure that the neurostimulator is placed no deeper than 2.5 cm (1 in) below the skin and is parallel to the skin. If the neurostimulator is too deep or is not parallel to the skin, telemetry may be unsuccessful.
- Do not coil excess length in front of the neurostimulator. Wrap excess length around the perimeter of the neurostimulator (Figure 5) to help minimize potential damage during neurostimulator replacement surgery, help minimize potential kinking, and minimize interference with telemetry during programming.


Figure 5. Wrap excess length around perimeter of the Model 3058 Neurostimulator.

2. Check the system integrity before securing the neurostimulator in place.
 - a. Use the clinician programmer and the product literature packaged with it to confirm the integrity of the connected system.

 **Warning:** To use the nonsterile programming system components in a sterile field, place a sterile barrier between the patient and programming system components to prevent infection. Do not sterilize any part of the clinician programmer. Sterilization may damage the programmer.

Note: The neurostimulator should be in the pocket during system interrogation to ensure proper readings.

- b. Program the stimulation parameters you have selected for your patient according to the product literature packaged with the programmer.
 - c. Check the battery status; if applicable, check the electrode impedances to rule out a short or open circuit.
3. Use the two suture holes in the connector block to secure the neurostimulator to the muscle fascia with nonabsorbable silk.

Completing the implant procedure

1. Close and dress all incisions.
2. Ensure that a patient control device and a patient ID card are given to the patient.

Caution: Because the patient control device is the patient's only means to adjust or turn the neurostimulator on or off, the patient must carry a patient control device at all times.

3. Schedule regular patient follow-up appointments to monitor the condition of the neurostimulator and to confirm that the programmed parameter values are appropriate.

Replacing the InterStim Neurostimulator

1. Open the implant site using normal surgical procedure and carefully remove the neurostimulator from the subcutaneous pocket.

Caution: If any component shows signs of damage, pitting, or corrosion, do not allow it to remain implanted. Damaged devices may fail to operate properly which could require additional surgery.

2. Clean the neurostimulator connector block and lead with sterile water; wipe lead dry with sterile gauze.
3. Insert a torque wrench through each prepierced hole in the rubber sealing grommet and loosen the setscrew(s) by turning it counterclockwise.
4. Gently retract the lead or extension from the neurostimulator connector block.
5. Clean and dry the lead, which must be free of fluids or tissue.

Note: Disconnect the extension (if applicable) from the lead. Clean and dry the proximal end of the lead. Take care not to move the lead when extension is disconnected.

6. Set aside the explanted components for return to Medtronic.
7. Follow the instructions for implanting the replacement neurostimulator according to the product literature packaged with the device.
8. Return explanted devices to Medtronic according to product literature packaged with those devices.

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
- **System overview and compatibility insert - Medtronic sacral neuromodulation therapy**

Ashley Luloff
Pr. Regulatory Affairs Specialist
Medtronic Inc.

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Asia-Pacific
Medtronic International Ltd.
50 Pines Road, #04-01,
Singapore Business City,
Singapore 117384,
Singapore
Tel: +65-6870-5510

Medtronic and Medtronic logo are
trademarks of Medtronic. All other marks
are trademarks of their respective companies.

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
710 Medtronic Pkwy
Minneapolis, MN 55432
USA
www.medtronic.com

Declaration of Authenticity

Medtronic Inc. located at 710 Medtronic Parkway, Minneapolis, MN 55432, USA hereby declares that the following document is true and accurate copy of the original kept on file at Medtronic Inc.


- **Clinician application for sacral neuromodulation therapy - Clinician Programming Guide for InterStim Systems**

Ashley Kuloff
Pr. Regulatory Affairs Specialist
Medtronic Inc.

Date: 09 January 2023.

Notary

Kimberly J. Horsch
January 10,

CLINICIAN APPLICATION FOR SACRAL NEUROMODULATION THERAPY

Clinician Programming Guide for
InterStim™ Systems

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Medtronic

SYMBOLS

Explanation of symbols in this manual

Refer to the appropriate product for symbols that apply.


Manufacturer

Conformité Européenne (European Conformity).
This symbol means that the device fully
complies with applicable European Union Acts.

Authorized representative in the European
community

For USA audiences only

Magnetic Resonance (MR) Conditional

Importer

Medtronic and the Medtronic logo are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company.

Android™ is a trademark of Google LLC.

The handset, as configured, does not support voice communication and therefore cannot be used to make any calls (including emergency calls).

Information available for the system:

The information for prescribers manual provides information about contraindications, warnings, precautions, adverse events, sterilization, patient selection, individualization of treatment, and component disposal.

For customers in Japan, the appropriate package insert provides information about indications, safety, contraindications, warnings, precautions, and adverse events.

The indications sheet provides information about indications and related information. For customers in Japan, the appropriate package insert provides information about indications.

The system eligibility, battery longevity, specifications manual provides information about neurostimulator selection, battery longevity calculations, and battery characteristics.

The system overview and compatibility insert provides information about component compatibility.

MRI guidelines provide information about any MRI conditions and MRI-specific contraindications, warnings, and precautions for MRI scans with the neurostimulation system.

Product manuals, such as programming guides, recharging guides, and implant manuals provide device descriptions, package contents, device specifications, product-specific warnings and precautions, and instructions for use.

! USA The clinical summary provides information about the clinical study results for the neurostimulation system.

Refer to the literature provided by the clinician handset manufacturer for information regarding wireless use.

Refer to the Communicator technical manual for warnings, precautions, device description, package contents, device specifications, instructions for use, and maintenance information.

TABLE OF CONTENTS

GENERAL WARNINGS AND CAUTIONS..	10
OVERVIEW	11
Intended Use	11
About This Guide	11
QUICK TIPS	12
Using the Workflow Navigator	12
Indicators and Icons	12
GETTING STARTED	14
Accessing the Clinician App	14
Creating a Password.....	15
Menu Screen	15
Demo Mode	16
Using the Communicator to Pair the Handset and Neurostimulator	16
Setting Implant Date.....	17
WORKFLOWS.....	19
Starting a Workflow.....	19
Configuring the Implanted Neurostimulator (Configured Implant Workflow)	20
Setting up implant configuration	20
Entering patient information	20
Checking device and system status.....	21
Therapy Parameter Settings	26
Maximum settings (Model 97810 InterStim Micro rechargeable neurostimulator only)	33
Optimizing device settings.....	33

Selecting Cycling feature and time.....	34
Performing an Integrity Measurement (Check Impedance Workflow)	37
Electrode polarity	38
Checking Neurostimulator MRI Eligibility (View MRI Status Workflow)	39
MRI - PREPARING FOR AN MRI SCAN	41
WORKING WITH PROGRAMS	46
Standard Programs	46
Custom Programs	46
Enabling programs for patient use.....	50
Reviewing a Therapy Programming Session ...	50
Ending a programming session.....	51
Using the Stop Therapy feature.....	51
DIAGNOSTICS	52
Usage Graph	52
To view the usage graph:	54
Event Log	54
Usage report	56
Recharge sessions (Model 97810 only)	56
REPORTS	57
Accessing a Downloaded Report	57
ADDITIONAL FEATURES	59
About Screen	59
Managing Devices	59
Model A901 Communication Manager Application	

60

Model A902 Patient Data Service Application .. 60

Software Information 60

Data security 60

Network connectivity 60

Installing application updates 61

Installing communicator updates 61

Returning the handset 62

TROUBLESHOOTING 63

Maximum Settings Reached (Model 97810 only) 73

Considerations for troubleshooting 73

Clinician Reset and Recharge (Model 97810 only)

73

Performing Clinician Reset and Recharge..... 73

GENERAL WARNINGS AND CAUTIONS

GENERAL WARNINGS AND CAUTIONS

WARNINGS:

Wound contact — Do not use the communicator on an unhealed wound. The communicator is not sterile, and contact with the wound may cause an infection. Keep a sterile bandage or barrier between the wound and device.

Sterile field — When using the communicator in a sterile field, place the communicator in a sterile bag. The communicator is not sterile and cannot be sterilized.

CAUTIONS:

Telemetry signal disruption from EMI — Do not attempt telemetry near equipment that may generate electromagnetic interference (EMI). EMI may cause a disruption in clinician programmer function.

Unresponsive app or handset — Use the recommended troubleshooting steps if the clinician app or handset becomes unresponsive. An unresponsive app or handset may result in an inability to complete the implant procedure or other programming sessions.

Note: Other **Warnings** and **Cautions** are found throughout this manual.

OVERVIEW

The Medtronic Model A51200/A510 clinician applications (apps) are intended for use with the HH90 Handset and TM90 Communicator to program, adjust, and troubleshoot the Medtronic Model 97810 InterStim Micro rechargeable neurostimulator and Models 3023 and 3058 InterStim neurostimulators for sacral neuromodulation therapy. The clinician uses the appropriate clinician app to program settings for the patient.

Intended Use

The clinician app is intended to program neurostimulators for sacral neuromodulation therapy.

About This Guide

This programming guide contains information a clinician needs to set up, use, and troubleshoot the clinician app on the handset. This guide provides instructions on how to use the clinician app to program an implanted neurostimulator. This guide is intended to be used for InterStim system implantable neurostimulators.

- Before using the clinician app for the first time, refer to the information provided with the handset for initial handset setup.
- All patient-related instructions on using the patient app on the handset are included in the A52200/A520 Patient Therapy App Patient User Guide for InterStim systems.
- Patient-related instructions on using the communicator are included in the TM90 Communicator Instructions for Use provided with the device.
- General instructions on recharging the Model 97810 InterStim Micro neurostimulator battery are included in the WR9220 Recharger Patient User Guide provided with the recharger.

QUICK TIPS

See the following information for helpful tips to aid you in using the clinician app, and to help you better understand the contents of this programming guide.

Using the Workflow Navigator

The workflow navigator is displayed at the top of every screen when you are in a workflow and represents the screens/ tasks that can be completed on each screen. A title that is highlighted and underlined in the workflow navigator indicates your position in that workflow. For more information on workflows, see "WORKFLOWS" on page 19.

You can use the workflow navigator in the following ways:

- Swipe to the left or right to view previous or subsequent screens within a workflow.
- Tap on a title to jump to a different screen.

A title that is highlighted and underlined indicates your position within a workflow.


Figure 1. Workflow navigator

Indicators and Icons

The following is a list of the icons found in the clinician app along with their associated descriptions.

Notes:

- The following is not an all-inclusive list of all the icons found in the clinician app, but rather a list of the more

important icons or those that appear most frequently within the app.

- If an icon or button appears gray within the clinician app, that option is not available.

Icon	Icon description
	Tap for additional screen options.
	Indicates that the neurostimulator device has been previously paired.
	Tap to access the About screen, Manage Devices screen, or end a programming session. This icon appears only when in a workflow or on the Home screen.
DEMO:	Indicates that the clinician app is operating in Demo mode.
	Tap to check battery level.
	Tap to check impedance.
	Tap to exit the current screen and go back to the one you were previously on.
	Indicates the program is enabled for patient use, and is visible to the patient.
	Indicates the program is not enabled for patient use, and is not visible to the patient.
	Tap to end the current programming session and exit the app.

GETTING STARTED

Before using the clinician app to program the neurostimulator, see the following sections for basic information about:

- Accessing the clinician app
- Creating a password
- Using the menu screen
- Operating in Demo mode
- Pairing to the implanted neurostimulator

Accessing the Clinician App

You can access the clinician app by locating the appropriate icon on the main screen of the handset (Figure 2). To open and begin using the clinician app, tap the Micro Clinician icon for the Model 97810 InterStim Micro rechargeable neurostimulator, or the Clinician icon for the Models 3023 and 3058 InterStim neurostimulators. Throughout this manual, “clinician app” refers to both apps, unless stated otherwise.

Figure 2. A51200 Micro Clinician app for the InterStim Micro rechargeable neurostimulator (left) and A510 Clinician app for Models 3023 and 3058 InterStim neurostimulators (right)

Note: The first time you open the app, it will request permission for data storage and session overlay. Tap **OK** to continue. Then tap **ALLOW** on the next message.

Creating a Password

The clinician app requires a password to use the app. A password is required to be entered under the following conditions:

- Upon launch of the app
- After 30 minutes of inactivity

Note: When you launch the clinician app for the first time, you will be asked to enter the default password. Contact Medtronic support for the default password. You will need to obtain the default password before creating your own unique password.

To create a password:

1. Enter the Medtronic default password in the appropriate field. Contact Medtronic support to obtain the default password if you haven't already.
2. Tap **OK** to continue to the next screen.
3. Enter your unique password in the appropriate field.

Note: The password is case-sensitive and must contain at least six characters.

4. Confirm your new password by entering it a second time in the appropriate field and tap **CREATE** when finished.

Note: For information on resetting a forgotten password, see "TROUBLESHOOTING" on page 63.

5. Select your region using the options provided on the handset.
6. Tap **NEXT** to continue on to the Menu screen.

Menu Screen

Upon entering your password, you are taken to the Menu screen. On the Menu screen, you can access Demo mode, connect the communicator and neurostimulator to the handset, or access patient reports.

Demo Mode

Demo mode allows you to simulate using the clinician app without being paired to an actual neurostimulator device. DEMO is displayed next to the screen title when the clinician app is operating in Demo mode.

Using the Communicator to Pair the Handset and Neurostimulator

All programming functionality takes place within the Connect option on the Menu screen. Tapping this option allows you to begin the process of using the communicator to pair the implanted neurostimulator to the clinician app on the handset.

To pair the handset and implanted neurostimulator using the communicator:

1. Press  on the communicator to turn the communicator on.
2. Navigate to the Menu screen, and tap **CONNECT**.
3. Place the communicator near the handset.
4. Tap  if you see this icon. The communicator will attempt to connect to the handset. The blue LED indicator on the communicator will continuously blink to indicate it is on and in discovery mode.

Note: If you do not see  proceed to step 5.

5. If you have previously paired the communicator to the handset, the communicator will automatically attempt to connect to the handset. If the blue LED indicator on the communicator is solid and no longer blinking, the communicator is now connected to the handset, skip to step 8.

Note: If you are having trouble connecting the communicator to the handset, see "TROUBLESHOOTING" starting on page 63.

6. If the communicator has not been previously paired to the handset, tap **SWITCH COMMUNICATOR** on the

handset.

7. Select the appropriate communicator device according to the serial number shown on the Connection screen, and tap **CONTINUE**. The blue LED indicator on the communicator will be solid and no longer blinking. This signifies the communicator is now connected to the handset.

Note: The serial number, e.g. NPA-xxxxxxx, is located on the back of the communicator, and should match what is listed on the Connection screen.

 **CAUTION:** Ensure you are connecting to the proper communicator by confirming the serial number on the device. Failure to do so could extend the implant procedure or programming session.

8. Ensuring the communicator is on, place the communicator over the implanted neurostimulator site and tap **FIND DEVICE** on the handset.
9. Select the appropriate neurostimulator device according to the serial number and tap **CONTINUE**. If the neurostimulator serial number is not listed on the Connection screen, tap **RETRY**.

Notes:

-  displayed next to the serial number of the neurostimulator device indicates which device has been most recently paired to the handset.
- Connecting to a new INS may take some time. If the INS has not finished connecting after a few minutes, press **Cancel** and try again.
- If you are having problems pairing to the neurostimulator, refer to "TROUBLESHOOTING" on page 63.

Setting Implant Date

Set implant date in the field shown on the screen in Figure 3.
To view implant date and Elective Replacement Indicator

CONTINUED

(ERI) date, navigate to the **About** screen. (Model 97810 InterStim Micro rechargeable neurostimulator only.)


Figure 3. Set Implant Date screen (Model 97810 InterStim Micro rechargeable neurostimulator only)

WORKFLOWS

Workflows are used to assist you in navigating programming sessions. Workflows for programming an implanted neurostimulator include:

- Configure implant
- Check impedance
- View MRI status

Starting a Workflow

Starting a workflow allows you to complete tasks for a programming session. You can access workflows on the Home screen.

To start a workflow:

1. Tap the name of a workflow to see an expanded view and description.
2. Tap **START** to proceed with the desired workflow.

Configuring the Implanted Neurostimulator (Configure Implant Workflow)

The Configure Implant workflow allows you to configure the neurostimulator for a patient. The main features of the Configure Implant workflow include:

- Manage implant configuration
- Manage patient information
- Configure therapy parameters
- Check therapy and/or electrode impedance
- Check battery status
- Run a battery longevity estimate (Model 3023 and 3058 neurostimulators only)
- Enable programs for patient use
- Create custom programs
- View a diagnostic report
- View the summary screen

Setting up implant configuration

The Configure Implant screen allows you to enter information required for determining the patient's MRI eligibility, such as lead model, device implant location, lead tip location, and the presence of abandoned product from removed systems. Even partial systems can have MRI interactions such as RF heating.

 **WARNING:** Ensure the implant configuration is up-to-date and accurate. If component information is entered incorrectly, the patient could be unnecessarily restricted from having an MRI scan or could be allowed to have an inappropriate MRI scan, which could cause serious or permanent patient injury.

Entering patient information

The Patient screen allows you to enter optional patient information (name, patient ID, gender, and date of birth).

Checking device and system status

The Status screen provides information about the neurostimulator and allows you to check the battery level of the neurostimulator and system impedance. Notice the battery and impedance rows; tap the battery and impedance icons for details. Table 1 displays the device status icons for impedance check and their associated description and meaning.

Table 1. Neurostimulator impedance check icons and descriptions

Icon	Description	Meaning
	Orange circle with exclamation point	Status is suspect. Requires further investigation.
	Green circle with check mark	Status is good. No further investigation needed.

Note: A question mark appears in the impedance row if an impedance check has not been performed. Electrodes with impedance measurements outside the nominal impedance range will display an orange icon.

Checking battery level

Tables 2 and 3 show the battery display of InterStim neurostimulators.

To check neurostimulator battery level for Model 97810 InterStim Micro rechargeable neurostimulator:

1. Navigate to the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Status** screen and tap  next to the battery level to see the current battery levels of the neurostimulator, communicator, and handset.
4. Tap **DONE** to exit the Battery screen when finished.

Table 2. Battery descriptions (InterStim Micro rechargeable neurostimulator only)

Approximate Battery Capacity	Model 97810 InterStim Micro rechargeable neurostimulator
Status: 30%-100% 	Battery Level: OK
Status: 20% 	Battery Level: LOW Recharge needed.
Status: 10% 	Battery Level: CRITICAL Recharge needed. Therapy may stop and the neurostimulator will need to be recharged soon. No access to the clinical application after you have received this notification.

Notes:

- Recharge the neurostimulator prior to implant if needed.
- Instruct patient to check rechargeable neurostimulator battery regularly to avoid loss of therapy.

To check neurostimulator battery level for Model 3023 and Model 3058 neurostimulators:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Status** screen to see the current battery level of the neurostimulator.
4. Tap **DONE** to exit the Battery screen when finished.

Table 3. Battery status descriptions (Models 3023 and 3058 neurostimulators only)

Approximate Battery Capacity	Model 3023 Neurostimulator	Model 3058 Neurostimulator
More than 15%	 Status: Battery Level: OK Longevity: Estimated number of months of battery longevity (based on new battery).	 Status: Battery Level: OK Longevity: Estimated number of months of battery longevity (based on estimated remaining capacity until LOW is displayed).
5-15%	 Status: Battery Level: LOW Battery is approaching End of Service, see below.	 Status: Battery Level: LOW Battery is approaching End of Service, see below.
Less than 5%	Notification: End of Service (EOS). Therapy will be stopped and the neurostimulator needs to be replaced. No access to the clinician application beyond this notification.	Notification: End of Service (EOS). Therapy will be stopped and the neurostimulator needs to be replaced. No access to the clinician application after you have received this notification.

Notes:

- Do not implant the neurostimulator if the battery status displays LOW or EOS prior to implant.
- If the neurostimulator battery level is below 15%, a low battery notification displays on the Home screen. The notification persists until the neurostimulator is replaced; however, you have the ability to dismiss the notification at any time.

Longevity estimate (for Models 3023 and 3058 InterStim neurostimulators only)

The longevity estimate provides therapy impedance, neurostimulator battery longevity estimate, and program settings for the active program. A longevity estimate for the implanted neurostimulator can be obtained from the Battery screen.

To obtain a neurostimulator battery longevity estimate:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the Status screen, and tap .
4. Tap **THERAPY MEASUREMENT**.
Note: Therapy must be turned on.
5. Tap **DONE** to return to the Battery screen.
6. Tap **DONE** again to return to the Status screen.

Notes:

- The longevity measurements are conservative estimates and are not intended to be device replacement indicators. Use LOW and EOS alerts for decisions on device replacement.
- If Cycling or SoftStart/Stop™ is enabled, longevity measurements may be inaccurate (Models 3023 and 3058 only). See "Selecting Cycling feature and time" on page 34 as well as the System Eligibility, Battery Longevity, Specifications Reference Manual for more information.
- Displayed battery longevity is an estimate, based on estimated battery capacity (Model 3058 only) and active program settings. Actual longevity will depend on any changes made to program settings and patient use throughout the service life of the INS. Also refer to the System Eligibility, Battery Longevity, Specifications Reference Manual for battery longevity calculations.

For the Model 3058 neurostimulator:

- The Model 3058 neurostimulator should be allowed to stabilize at body temperature for at least one week before conducting Therapy Measurements. The longevity measurements provided prior to this period or prior to implant may be lower than expected.
- If the Model 3058 neurostimulator has experienced a Power on Reset (POR) the estimated number of months of battery longevity will be calculated as if the battery was new.
- If the estimated battery capacity is greater than 15% of the original capacity, the following additional information is displayed:
 - Number of months based on remaining estimated battery capacity and program settings.

Impedance measurements

You can check impedance on the Status screen of the Configure Implant workflow, or within the Check Impedance workflow.

Note: Electrodes with impedance measurements outside the nominal impedance range will display an orange icon. See "Table 1. Neurostimulator impedance check icons and descriptions" on page 21 for details.

To check impedance from the Configure Implant workflow:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **Start**.
3. Navigate to the **Status** screen.
4. Tap .
5. Tap  again to initiate an impedance check.
6. Tap **DONE** when finished.

Notes:

- Tapping  on the status screen will take you to the Impedance screen where you can view the current impedance values of the respective electrode pairs and neurostimulator case.
- You can tap the status icon on any electrode to view the current impedance values.
- For Models 3023 and 3058 neurostimulators, you have the option to adjust therapy amplitude and/or pulse width temporarily while conducting an impedance check on the Impedance screen. Tap  in the corner of the screen to temporarily adjust amplitude and adjust pulse width parameters.

Therapy Parameter Settings

You can adjust therapy parameter settings as needed for a program in order to provide optimal therapy. You can access therapy parameter settings by navigating to the **Therapy** screen in the Configure Implant workflow.

During a programming session, the clinician app modifies the specified therapy parameters to program the neurostimulator. Table 4 lists the available therapy parameter settings along with their descriptions.

Table 4. Therapy parameter descriptions

Therapy parameter setting	Description
Amplitude	Amplitude is the intensity or strength of the stimulation measured in milliamps (mA) or in volts (V). By increasing amplitude, you are increasing the intensity of the stimulation and by decreasing it, you are reducing the intensity. The goal is to use the lowest effective amplitude that will provide optimal patient symptom relief, minimize patient discomfort, and maintain neurostimulator battery life to the best possible extent.

Therapy parameter setting	Description
	A typical setting for amplitude is based on patient comfort, which is determined during a programming session. Available amplitude ranges are: - For the Model 97810 InterStim Micro rechargeable neurostimulator: 0-12.5 mA - For the Model 3023 neurostimulator: 0-10.5 V - For the Model 3058 neurostimulator: 0-8.5 V
Pulse Width	Pulse width is the time or duration of the stimulation pulse measured in microseconds (µsec). Increasing pulse width increases pulse duration and decreasing pulse width decreases pulse duration. For example, when a patient feels the stimulation too intensely in one body location, increasing the pulse width spreads the stimulation and makes it less intense in that location. Default setting for pulse width is 210 µsec.
Rate	Rate is the number of times per second a pulse is delivered, measured in pulses per second (pps) or Hertz (Hz). Increasing the rate feels more like a "flutter" or "vibration" and decreasing the rate gives more of a "tapping" or "thumping" sensation. You can use patients' preferences for the sensation they are most comfortable with, to guide you in selecting an appropriate rate. Default setting for rate is 14 pulses per second.
Optional settings: The following are optional therapy parameters that can be modified if optimal therapy is not provided through default standard programs.	
Amplitude limit	Selecting an Amplitude Limit sets the boundary for how high the amplitude value can be raised above the programmed therapy amplitude value by the patient. After a programming session, patients can control therapy amplitude using the patient app, but within the programmed Amplitude Limit set by the clinician. When Amplitude Limit is not defined, patients can increase therapy amplitude to the maximum therapy amplitude value allowed.

Therapy parameter setting	Description
Cycling	Cycling is a feature used to control how often therapy is running. When the Cycling feature is off, therapy is continuously running. When the Cycling function is turned on, it automatically turns the neurostimulator on and off at clinician-determined intervals (from 2 seconds to 24 hours for Model 97810 and from 0.1 seconds to 24 hours for Models 3023 and 3058). When Cycling On time is over, therapy turns off and continues to stay off until Cycling Off time ends. When Cycling Off time is over, therapy turns on again and continues to stay on until Cycling On time ends. This process repeats itself until the Cycling feature is disabled or therapy is stopped. By default, the Cycling feature is turned off. Note: Due to a carryover effect, the patient may continue to experience symptom improvement during the Cycling Off time.
SoftStart/Stop	This feature is intended to increase patient comfort by providing a gentle or "soft" start as stimulation begins and reduces the risk that the patient will be startled by the start of a stimulation cycle. Using the SoftStart/Stop feature helps avoid any unpleasant sensation at the onset of stimulation. Programmable ramp times for this feature (up and down) are 1, 2, 4, 8, 15, or 30 seconds. By default, SoftStart/Stop is turned On with a ramp time of 4 seconds.
Amplitude resolution (Models 3023 and 3058 only)	Amplitude resolution is a feature that allows you to change the increments by which amplitude is increased or decreased, and allows you to ramp up or down to the appropriate amplitude setting. The amplitude resolution is provided through two options: • Default • Fine Changing amplitude from default to fine resolution is beneficial if a patient is sensitive to the intensity, measured in Volts (V), of the amplitude being provided. If the default resolution is selected, amplitude increases/decreases by increments of 0.1V. If the fine resolution is selected, amplitude increases/decreases by increments of 0.05V. Range for both Models 3023 and 3058: 0-6.35V.

Therapy Parameter Settings Instructions For Use

To turn on therapy by setting the amplitude:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Use the **arrows** to set the amplitude. Hold the Up or Down arrows to adjust the amplitude by increments of 0.5 mA or V.

Note: If the patient is experiencing discomfort, therapy can be stopped at any time. Tap  to stop therapy.

5. Tap **DONE** when finished.

Indicates that the program is visible to the patient.

Displays the electrode configuration for the program.


Figure 4. Therapy screen

Note: Model 97810 InterStim Micro rechargeable neurostimulator shows amplitude as mA instead of V.

To set pulse width:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Select **Pulse Width** from the list of options.
6. Use the **arrows** to set the pulse width.
 **CAUTION:** Consider default settings for rate and pulse width when configuring therapy settings. Unsuitable stimulation could result in patient discomfort.
7. Tap **UPDATE**.
8. Tap **DONE** when finished.
9. Tap **DONE** again to return to the Therapy screen.

To set rate:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Select **Rate** from the list of options.
6. Use the **arrows** to set the rate.
 **CAUTION:** Consider default settings for rate and pulse width when configuring therapy settings. Unsuitable stimulation could result in patient discomfort.
7. Tap **UPDATE** when finished.
8. Tap **DONE** when finished.
9. Tap **DONE** again to return to the Therapy screen.

To view electrode settings:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Select **Electrodes** from the list of options. The Electrode screen displays the current electrode settings for the selected program.
6. Tap **DONE** when finished.

Notes:

- Electrode settings for programs 1-7 are set up by default and cannot be changed.
- You can also view electrode settings for programs on the Therapy screen (Figure 4).

To set amplitude limit:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Select **Amplitude Limit** from the list of options.
6. Tap the **On/Off** switch to enable amplitude limit.

Note: Tap the **On/Off** switch again to disable.

7. Use the **arrows** to set the amplitude limit.
8. Tap **UPDATE** when finished.
9. Tap **DONE** when finished.
10. Tap **DONE** again to return to the Therapy screen.

Notes:

- When the On/Off switch is turned on, the amplitude limit will default to the current amplitude setting.
- The amplitude limit cannot be set at a value lower than the current amplitude setting for an active program.

Maximum settings (Model 97810 InterStim Micro rechargeable neurostimulator only)

The system cannot use certain combinations of amplitude, pulse width, and rate when the battery is low, because the neurostimulator has reached the limit of its therapy range. If this happens, a Maximum Settings Reached alert will appear, and an orange icon will appear to the left of the parameter value, indicating that stimulation is not being provided at the level that is shown on screen.

To resolve this condition, perform one of the following actions:

- Slowly reduce the stimulation parameter until the icon disappears.
- Reduce other parameter values to set a desired value higher. For example, if high amplitude is desired, reduce the rate or the pulse width.
- Check the battery level. If the battery level is low, consider recharging the neurostimulator

For additional information on troubleshooting these types of errors, refer to "Maximum Settings Reached troubleshooting" on page 73.

Optimizing device settings

The clinician app offers features that you can optionally use to increase patient comfort and when certain parameters

are used, extend neurostimulator battery life. These features include:

- Cycling
- SoftStart/Stop
- Amplitude resolution (for Models 3023 and 3058 neurostimulators only)

See "Therapy Parameter Settings" on page 26 for more information on setting additional therapy parameters.

Selecting Cycling feature and time

The Cycling feature automatically turns the neurostimulator On and Off at clinician-determined intervals as follows:

- From 2 seconds to 24 hours for Model 97810
- From 0.1 seconds to 24 hours for Models 3023 and 3058

Cycling, by default, is disabled.

For Models 3023 and 3058 neurostimulators only: Anytime Cycling is enabled, the battery longevity estimate may not be accurate. If Cycling is being used to improve device battery longevity, the information below should be considered:

- For amplitude settings less than or equal to 1 Volt, Cycling may result in reduced device battery longevity when compared to continuous mode.
- For amplitude settings greater than 1 Volt, the following settings may improve device battery longevity when compared to programming in continuous mode:
 - Set Cycling intervals to greater than or equal to 2 seconds On and greater than or equal to 2 seconds Off (without SoftStart/Stop enabled).

Notes:

- When using both SoftStart/Stop and Cycling features, only SoftStart/Stop times that are less than or equal to the Cycling On and Cycling Off times are available for programming.

 **CAUTION:** In Models 3023 and 3058 neurostimulators, using some Cycling or SoftStart/Stop therapy settings may cause premature battery depletion, resulting in loss of therapy.

To enable SoftStart/Stop:

1. Navigate to the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Tap **SoftStart/Stop** from the list of options.
6. Tap the **On/Off** switch to enable SoftStart/Stop.

Note: To disable, tap the **On/Off** switch again.

7. Use the **arrows** to configure the SoftStart/Stop time for the implanted neurostimulator.

Note: For Models 3023 and 3058 neurostimulators only, the SoftStart/Stop time cannot be more than the Cycling On or Cycling Off time, and/or the Cycling On time and the Cycling Off time cannot be less than the determined SoftStart/Stop time.

8. Tap **UPDATE**.

Note: Certain Cycling or SoftStart/Stop configurations may reduce device battery longevity. (Models 3023 and 3058 neurostimulators only)

9. Tap **DONE** when finished.
10. Tap **DONE** again to return to the Therapy screen.

To enable Cycling:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program

you would like to configure (Figure 4).

4. Tap .
5. Tap **Cycling** from the list of options.
6. Tap the **On/Off switch** to enable Cycling.

Note: Tap the **On/Off switch** again to disable.

7. Use the **arrows** to set the Cycle On time and Cycle Off time for the implanted neurostimulator. The length of Cycle On time and Cycle Off time is measured in seconds.

Note: Certain Cycling or SoftStart/Stop configurations may reduce device battery longevity. (Models 3023 and 3058 neurostimulators only)

8. Tap **UPDATE**.
9. Tap **DONE** when finished.
10. Tap **DONE** again to return to the Therapy screen.

To change amplitude resolution (Models 3023 and 3058 neurostimulators only)

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure (Figure 4).
4. Tap .
5. Tap **Resolution** from the list of options.
6. Select the desired amplitude resolution (default or fine) from the available options.
7. If the amplitude resolution is at the desired setting, tap **DONE**.
8. If you have changed the amplitude resolution, once on the desired setting, tap **UPDATE**.
9. Tap **DONE** when finished.
10. Tap **DONE** again to return to the Therapy screen.

WORKFLOWS

Performing an Integrity Measurement (Check Impedance Workflow)

An impedance measurement can be run to check system integrity. Run an impedance check at the beginning of every programming session.

To check impedance from the Check Impedance workflow:

1. Select the **CHECK IMPEDANCE** workflow.
2. Tap **START**.
3. Tap  to check electrode impedance (Figure 5).
4. Tap **DONE** when finished.

Notes:

- Tapping  on the status screen will take you to the Impedance screen where you can view the previous impedance measurement of the electrodes.
- You can tap the status icon on any electrode to view the individual impedance values.
- You have the option to adjust therapy amplitude and/or pulse width temporarily while conducting an impedance check on the Impedance screen. Tap  in the corner of the screen to access temporary adjust amplitude and adjust pulse width options. (Models 3023 and 3058 neurostimulators only)

CONFIDENTIAL

Tap the numbered circle icons to view impedance measurement values for the corresponding electrode.


Figure 5. Impedance screen

Note: Only Models 3023 and 3058 neurostimulators show  in the corner of the Impedance screen.

Electrode polarity

The stimulation pulse is delivered from the neurostimulator to the nerve through the electrodes on the leads. For the stimulation pulse to reach the nerve, you select the

electrodes on the lead that best provide the stimulation and assign a negative or a positive polarity. At least one electrode (or the neurostimulator case) must be designated as positive and at least one electrode must be designated as negative. The negative electrode is called the active electrode; a pulse flows from the active or negative electrode to the positive. Thus, changing an electrode to active changes the location of the stimulation pattern.

Configuring electrodes depends on how the leads and electrodes are placed in relation to the nerve that is being stimulated. Approach configuration systematically, using a variety of electrode configurations and mapping sensory responses to each configuration. The four electrodes on a lead can be configured with or without the neurostimulator case:

- **Unipolar**—Using any combination of electrodes with at least one electrode as negative and the neurostimulator case as positive.
- **Bipolar** —Using any combination of electrodes with at least one electrode as positive, one negative, and the case off.

Notes:

- The neurostimulator case can only be configured as positive. When the case is selected, lead electrodes can only be selected as negative.
- Unipolar configuration depletes the neurostimulator battery at a higher rate. Use bipolar configuration whenever possible to extend battery life.

Checking Neurostimulator MRI Eligibility (View MRI Status Workflow)

The Model 97810 InterStim Micro rechargeable neurostimulator with SureScan™ MRI technology and Models 3023 and 3058 neurostimulators may be eligible for MRI scans. Follow the instructions beginning on page 41

to assess neurostimulator eligibility and to prepare the
neurostimulator for an MRI scan.

MRI - PREPARING FOR AN MRI SCAN

MRI - PREPARING FOR AN MRI SCAN

 **Magnetic Resonance (MR) Conditional:** Depending on what kind of neurostimulation system components the patient has implanted, the patient may be eligible for one of the following types of magnetic resonance imaging (MRI) scans:

- MRI scans of the head only (this is referred to as head-scan only eligible).
- MRI scans of any part of the body (this is referred to as full-body scan eligible).

Use the clinician app on the patient's handset to activate MRI mode, assess MRI eligibility, and prepare the patient's InterStim system for an MRI scan. For the MR conditions and MRI-specific warnings and precautions for conducting an MRI scan, refer to the MRI Guidelines for InterStim systems Manual.

 **WARNING:** Do not conduct an MRI scan on a patient before assessing MRI eligibility. Instruct the patient to inform the MRI clinician that the patient has an implanted device, and to activate MRI mode before the scan. Conducting an MRI scan on a patient who is ineligible for an MRI scan could result in device damage or serious or permanent patient injury.

 **CAUTION:** If the patient has a rechargeable neurostimulator battery, ensure the battery level of the neurostimulator is at least 30% before the patient undergoes an MRI scan. Patients are not able to receive an MRI scan if a rechargeable neurostimulator battery is not sufficiently charged.

 **CAUTION:** Deactivate MRI mode when the MRI scan is complete and the patient is outside the MRI scanner (magnet) room. Failure to deactivate MRI mode and return amplitude to the previous therapy level may result in a return of symptoms.

MRI - PREPARING FOR AN MRI SCAN

Instruct the patient prior to the MRI scan to:

- Bring the handset and communicator to the MRI appointment, however not to take any patient control device into the MRI scanner (magnet) room.
- Bring the patient identification (ID) card provided by Medtronic for the implanted neurostimulation system.
- Tell the MRI clinician that they have an implantable device.
- Ensure the battery status of the Model 97810 InterStim Micro rechargeable neurostimulator is at least 30% before undergoing an MRI scan.
- Activate MRI mode before MRI scan, outside the MRI scanner (magnet) room. MRI mode turns therapy off and prepares the system for MRI scanning. Because therapy is off, symptoms may return.

Notes:

-If the patient has an optional control magnet (Model 3023 neurostimulator only) disable the magnet switch before the MRI scan using the Model 8840 Clinician Programmer.

-Even though the patient's implanted system is eligible for certain MRI scans, the patient's scan-type eligibility depends on factors including:

- The area of the patient's body that needs to be scanned
- The ability of the MRI facility to meet all the required conditions for scanning

Viewing the patient's MRI scan-type eligibility and activating MRI mode:

1. From the Home screen, tap **VIEW MRI STATUS**.
2. From the Eligibility Identification screen, confirm that the information for the patient's implanted components (ie, neurostimulator and lead model) and their location are current and accurate. Tap **CONFIRM**.

Notes:

MRI - PREPARING FOR AN MRI SCAN

If the information for the implanted components and their location is not up-to-date and accurate, the resulting MRI scan-type eligibility may be incorrect.

If the handset screen says **Not Ready for MRI Scan** and **Not Eligible**, the system is not eligible for MRI.

3. Place the communicator over the device and tap **ACTIVATE** to activate MRI mode. Activating MRI mode turns therapy off.
4. Determine which one of the following messages appears on the MRI eligibility screen:
 - MRI Conditional Full Body Scan Eligible
 - MRI Conditional Head Scan Eligible with Transmit/Receive Head Coil
 - MRI eligibility cannot be determined


Note: If the handset screen says MRI eligibility cannot be determined, further assessment is necessary. Review reasons for result on screen.

5. See Table 5 for examples of the MRI eligibility messages. The screen confirms the patient's MRI scan-type eligibility and whether MRI mode is activated.
6. If you will not be present at the MRI facility or are not activating MRI mode prior to the patient's MRI scan, explain to the patient how to activate MRI mode using the handset at the MRI exam.

Note: You may also complete the MRI scan-type eligibility form provided in the MRI Guidelines for InterStim systems Manual and give the form to the patient so that they can take it to their MRI exam.

MRI - PREPARING FOR AN MRI SCAN

Table 5. SCAN ELIGIBILITY screen displays

Scan Eligibility	Screen image, text and description
MR Conditional Full Body Scan Eligible	MRI Mode is Activated  MR Conditional Full Body Scan Eligible
MR Conditional Head Scan Eligible with Transmit/Receive Head Coil	MRI Mode is Activated  MR Conditional Head Scan Eligible with Transmit/Receive Head Coil
Cannot be Determined	MRI Mode is Activated  MRI eligibility cannot be determined MRI mode is activated but scan-type eligibility cannot be determined. Review reasons for eligibility result on screen.
Device not eligible for MRI scan	Not Ready for MRI Scan  Not Eligible Model 3023 neurostimulator only: The neurostimulator may be ineligible for MRI, or the magnet switch inside the neurostimulator may be enabled. Before an MRI scan, the magnet switch for the Model 3023 neurostimulator must be disabled using the Model 8840 Clinician Programmer. Call Medtronic support for assistance.

Notes:

Refer to the MRI Guidelines for InterStim systems Manual for more information about determining MRI eligibility.

MRI - PREPARING FOR AN MRI SCAN

The "consult instruction for use" symbol () when shown with MRI scan eligibility means "consult the MRI guidelines for this neurostimulation system."

Deactivating MRI mode and turning therapy on after the patient's MRI scan:

1. Connect to the neurostimulator as usual.
2. From the Home screen, tap VIEW MRI STATUS.
3. Tap **DEACTIVATE** at the bottom of the screen.
4. Tap **DEACTIVATE** a second time when prompted.
5. Tap **YES** when prompted to return to previous therapy settings, if desired.

Note: If the Home screen is not shown after completing step 1, refer to "TROUBLESHOOTING" on page 63.

WORKING WITH PROGRAMS

Standard and custom programs are available within the clinician app. These features allow you to provide optimal therapy for the patient. Standard programs provide pre-defined electrode configurations. If default electrode settings from standard programs provide less than optimal therapy, custom programs allow you to manually define electrode settings.

Standard Programs

Programs marked 1-7 under Program Library on the Therapy screen of the Configure Implant workflow signify standard programs with pre-configured electrode settings. Only one program, either standard or custom, can be active at a time. The active program name is displayed at the top of the Therapy screen under **ACTIVE PROGRAM**.

Note: Upon switching from one active program to another, therapy amplitude automatically resets to zero.

Custom Programs

If the electrode settings for standard programs (1-7) do not provide optimal therapy results, you have the option to create custom programs using custom electrode settings. Programs A-D signify custom programs that you can create using custom electrode configuration settings. You can create up to four custom programs.

To create a custom program:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen.
4. Scroll to the bottom of the screen, and tap **+**.
5. Tap  to configure the electrodes for the custom program.
6. Set the electrode configuration by tapping in the **box**

for each electrode you want to configure. Tap once for a positive (+) electrode, tap again for a negative (-) electrode, and once more to clear (Figures 6 and 7).

Notes:

- Therapy must be off before modifying electrode configuration settings.
 - There must be at least one positive (+) electrode and one negative (-) electrode for a viable electrode configuration.
 - If the case is set to positive (+), the lead electrodes used must be negative (-).
7. Tap **UPDATE** after the electrode configurations have been set.
 8. Tap **DONE** when the update is complete.
 9. Adjust the parameter settings as desired. For step by step instructions, see "Therapy Parameter Settings Instructions For Use" on page 29 for more information.
 10. Tap **DONE** again to return to the Therapy screen.

Tap inside the individual numbered boxes to configure electrodes for custom programs.


Figure 6. Set electrode configurations for custom programs for the Model 97810 InterStim Micro rechargeable neurostimulator

2017-11-15 10:38 AM Page 49

Tap inside the individual numbered boxes to configure electrodes for custom programs.


Figure 7. Set electrode configurations for custom programs for Models 3023 and 3058 neurostimulators

Deleting a custom program

In the event that a custom program is no longer needed or used, it can be deleted at any time.

To delete a custom program:

1. Navigate to the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the custom program you would like to delete.
4. Scroll to the bottom of the screen, and tap **X**.

Note: You cannot delete a custom program if the custom program is active.

Enabling programs for patient use

In order for patients to use a program, the program needs to be enabled for patient use.

To enable a program for patient use:

1. Navigate to the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to enable.
4. On the bottom of the screen, tap the **Enable for patient use** switch once to make the program enabled for patient use.

Note: Tap the **Enable for patient use** switch again to disable the program for patient use.

5. Tap **DONE** when finished.  is displayed next to the program indicating that the program is visible to the patient.

Reviewing a Therapy Programming Session

The last screen in the Configure Implant workflow is the Summary screen. The Summary screen displays the active program setting and a list of patient enabled programs for a programming session.

Note: On the Summary screen, you can also view a session report. For more information, see "REPORTS" on page 57.

Ending a programming session

After completing the necessary programming tasks in a workflow, end the programming session.

 **CAUTION:** End programming session — Do not leave a programming session open in the clinician app once therapy has been configured. Not ending a programming session may result in patient access to the clinician app, resulting in the patient being able to configure inappropriate therapy.

To end a programming session:

1. Tap  in the corner of the screen.
2. Tap **End Session**.
3. Confirm that you would like to end the programming session by tapping **END SESSION**.

Note: Alternatively, you can also tap  on the **Home** screen to end a programming session.

Using the Stop Therapy feature

The Stop Therapy feature allows you to discontinue therapy at any time while in the Configure Implant workflow.

To stop therapy:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Therapy** screen, and tap the program you would like to configure.
4. Tap .
5. Tap **DONE** to return to the Therapy screen.

Note: For instructions on turning therapy back on, see "To turn on therapy by setting the amplitude:" on page 29.

DIAGNOSTICS

The clinician app allows you to view the following diagnostics to analyze a patient's therapy:


- Usage graph
- Event log
- Usage report
- Recharge session (Model 97810 InterStim Micro rechargeable neurostimulator only)

Usage Graph

The usage graph represents daily patient activity while using the therapy (Figure 8). The usage graph displays programs used, the amplitude setting and adjustments, as well as the duration of time each program was used. Each program is represented in the graph through the use of different colors and program numbers.

Time-stamp on the y-axis represents the duration of the program used.

Indicates the program used.

Amplitude values on the x-axis represent the amplitude value each program is set to.

Figure 8. Usage graph

Note: Model 97810 InterStim Micro rechargeable neurostimulator shows amplitude as mA instead of V.

To view the usage graph:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Diagnostics** screen, and tap .

Notes:

- Amplitude adjustments are only displayed for amplitude values that have been programmed for 30 minutes or more.
- The graph displays up to the last 30 days of patient data.

Event Log

The event log provides the following information (Figure 9):

- Date and time of therapy changes
- Amplitude adjustments (if programmed for 30+ minutes)
- Program changes
- Power-on-resets (POR)
- Change in therapy status (ON/OFF)

Figure 9. Event log

Note: Model 97810 InterStim Micro rechargeable neurostimulator shows amplitude as mA instead of V.

To view the Event log:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Diagnostics** screen.
4. Tap 

Usage report

The usage report displays data from both the usage graph and event log.

To generate a usage report:

1. Select the **CONFIGURE IMPLANT** workflow.
2. Tap **START**.
3. Navigate to the **Diagnostics** screen.
4. Tap **GENERATE USAGE REPORT**.

Note: You can access the generated report by navigating to **REPORTS** on the Menu screen.

Recharge sessions (Model 97810 only)

For the Model 97810 InterStim rechargeable neurostimulator, recharge sessions displays information about the patient's previous recharge sessions. This information includes:

- The date on which the recharge session occurred
- The duration of the recharge session
- The neurostimulator battery level at the start and at the end of the recharge session
- The average coupling strength

Note: The list of sessions can be sorted chronologically (either newest to oldest, or oldest to newest) with information for up to 35 of the most recent recharge sessions.

REPORTS

REPORTS

The clinician app allows you to generate reports that track data once the neurostimulator has been implanted and configured for chronic therapy. Access reports by tapping **REPORTS** on the Menu screen.

Table 6. Reports and descriptions

Report	Description
Session report	• Displays the patient's active program settings. • Generates automatically anytime the clinician ends a programming session.
Usage report	• Includes information from the event log and usage graph. • Generated from the Diagnostics screen.
Medtronic data report	• Generated from the Device tab on the About screen. • Used by Medtronic to troubleshoot the system.

To view, download, or delete reports:

1. Navigate to the **Menu** screen and tap **REPORTS**.
2. From the list of options, tap the type of report that you would like to view.
3. Tap the desired report to view a list of options.
 - a. Tap  to view a report.
 - b. Tap  to download a report.
 - c. Tap  to delete a report.

Accessing a Downloaded Report

If you need to download a report, see the following instructions on how to download the report and access it using your computer.

REPORTS

To download and access a report:

1. Navigate to the **Menu** screen and tap **REPORTS**.
2. From the list of options, tap the type of report that you would like to view.
3. Tap the desired report to view a list of options.
4. Tap  to download. A pop-up displays that the report has been saved to the Reports folder.
5. Plug the handset into a computer using the micro USB cable that came with the handset.
6. Navigate to the **Reports** folder from your computer drive. The Reports folder stores the reports that were downloaded from the clinician app.

Note: Generated reports will remain on the handset for 24 hours before being deleted automatically.

ADDITIONAL FEATURES

See the following sections for information on the additional features of the clinician app.

About Screen

You access the About screen by tapping  in the corner of the screen and selecting **About** from the list of options. The About screen provides the following information:

- Handset and app information such as model numbers, serial number, version numbers, etc.
- Neurostimulator device information such as model number and serial number.
- Communicator information such as firmware version number and serial number.
- General licensing information

Managing Devices

The Manage Devices screen allows you to view, add or remove a neurostimulator device. The Manage Devices screen can be accessed in the following ways:

1. From the Menu screen:
 - a. Tap .
 - b. Tap **Manage Devices** from the list of options.
2. After a device has already been connected:
 - a. Tap .
 - b. Tap **Manage Devices** from the list of options.
3. Tap **DONE** when finished.

Note: Tap  to remove a device, or  to add a new device. You will be redirected to either the Home screen or Menu screen depending on your selection.

Model A901 Communication Manager Application

The Model A901 Communication Manager app works with the Model A510 and A51200 Clinician Programming apps to manage telemetry communication between the handset and the neurostimulator. You do not need to open or use this app.

Model A902 Patient Data Service Application

The Model A902 Patient Data Service app works with the Model A510 and A51200 Clinician Programming apps to store patient and therapy device records that are used when generating a report.

Software Information

The following is important information regarding Medtronic therapy application software.

Data security

InterStim system clinician apps use and store data about the patient's health and implanted medical device. Data is protected by encryption within the application and encryption provided by the handset. The system does not protect data exported to another destination. Exported data should be handled in accordance with your facility's security policy for data handling and storage.

- Medtronic recommends that you always save exported data to the default reports location on the handset.

Network connectivity

Network connectivity is required for initial app registration and for installation of Medtronic app updates and communicator firmware updates. Network connectivity is not required for neurostimulator programming. To protect patient information,

Medtronic recommends you implement the following security measures:

- Secure the handset by disabling network connectivity during any programming session.
- Use a managed, trusted Wi-Fi connection when network connectivity is needed.

Note: Connecting the handset to a network that includes other equipment could result in unforeseen risks to patients, operators, or third parties. Changes to your network (such as adding, disconnecting, and upgrading equipment; upgrading or installing software; or changing network configurations) could also introduce additional risks.

- If your handset is lost or stolen, contact Medtronic right away using the contact information listed on the back cover of this manual.
- If you suspect a cybersecurity event has occurred, stop using the app (if possible), and contact your IT department or Medtronic support to document and respond to the suspected incident.

If you have any further questions related to cybersecurity, contact your IT department or Medtronic support.

Installing application updates

The clinician app along with other related apps are provided pre-installed on the Android® platform-based handset and cannot be installed on a personal cellular phone. Medtronic periodically updates the therapy app and will not install updates without notifying the user. Network connectivity is required to update the app. When notified that an app update is available, follow the instructions provided by Medtronic to install it.

Installing communicator updates

When you receive a notification to install a communicator firmware update, connect the handset to a network and follow the instructions.

Returning the handset

If you need to return the handset for disposal or replacement, contact Medtronic using the contact information listed on the back cover of this manual. Instructions will be provided for preparing the clinician handset for return.

TROUBLESHOOTING

The clinician app displays notifications you may encounter while using the app. If you are experiencing problems with the handset, check the app screen and follow the instructions to correct the issue. If you experience an issue not described in this section, contact your Medtronic support representative.

Table 7. Clinician app, handset, and communicator scenarios and solutions

Scenario	Explanations and possible solutions
I am unable to connect to the communicator	The communicator is not powered on. You will need to turn on the communicator to communicate with the implanted neurostimulator. Once the communicator is powered on, follow the connection steps to attempt to reconnect to an existing communicator. See, "Using the Communicator to Pair the Handset and Neurostimulator" on page 16 for instructions. Note: For information on features and functions of the communicator, refer to the TM90 Communicator Instructions For Use. The communicator is not charged or is in the process of charging. You will need to charge the communicator using the charger provided with the product package. Once charged, follow the connection steps to attempt to reconnect to an existing communicator. See, "Using the Communicator to Pair the Handset and Neurostimulator" on page 16 for instructions.

Scenario	Explanations and possible solutions
	The communicator is not paired to the handset. Follow the connection steps to attempt to reconnect to an existing communicator. See, "Using the Communicator to Pair the Handset and Neurostimulator" on page 16 for instructions. Note: The communicator must be paired to the handset first before attempting to connect to the neurostimulator. You received a replacement communicator. You will need to connect the new communicator to the handset and the implanted device. 1. Open the clinician app. 2. Make sure the communicator is on. Note: For information on features and functions of the communicator, refer to the TM90 Communicator Instructions For Use. 3. A screen will be displayed showing that the communicator is not found. Tap SWITCH COMMUNICATOR to connect. 4. Select the communicator you would like to pair to. The serial number of the communicator should match the serial number shown on the screen. 5. If the communicator you would like to connect to is not shown, tap RETRY. Follow the remaining instructions on the handset to finish pairing the device. 6. If the issue persists, contact Medtronic support.

TROUBLESHOOTING

Scenario	Explanations and possible solutions
I am unable to communicate with the implanted neurostimulator	The communicator is out of range. Reposition the communicator over the neurostimulator and tap RETRY. Make sure the communicator is flush with the skin. The rechargeable neurostimulator battery is depleted (Model 97810 InterStim Micro neurostimulator only). 1. Recharge the neurostimulator using the instructions provided with the recharging device. 2. If unsuccessful, follow the instructions for Clinician Reset and Recharge on page 73. 3. If the issue persists, contact Medtronic support.
Uncomfortable or intolerable therapy	The patient is experiencing side effects from the therapy (stimulation). 1. Turn the therapy (stimulation) down or off. See "Therapy Parameter Settings" on page 26 for instructions. 2. If the issue still persists, contact Medtronic support.
The handset has no power or has lost power	The handset battery is depleted. Recharge the battery using the charger. The handset is damaged or malfunctioning. Contact Medtronic support.
The handset will not charge	The charger is disconnected from the handset. Connect the charger to the handset.

Scenario	Explanations and possible solutions
	An incompatible charger is connected to the handset. Connect the appropriate charger to the handset. Charger is defective. Contact Medtronic support. The charger will need to be replaced. The handset is damaged or malfunctioning. Contact Medtronic support.
The patient is not receiving therapy	The therapy might be off. Turn on the therapy. The neurostimulator battery may be depleted. The neurostimulator will need to be recharged or replaced. A lead is damaged or has become disconnected. Check the impedance on the lead. If impedance is out of range, the lead may be damaged. If impedance is in range, select a different program. If the issue persists, contact Medtronic support.
I cannot find the clinician app on the handset	The clinician app may not have downloaded properly. Contact Medtronic support.

Scenario	Explanations and possible solutions
I have forgotten my password	The password needs to be reset. To reset the password 1. Tap the "Forgot Password" icon. 2. Take note of the six-character code. You will need to provide this to code to Medtronic support in order to complete the steps to reset your password. Note: This code is valid for 24 hours. After 24 hours, you will need to refresh the code. 3. Contact Medtronic support to complete the steps to reset your forgotten password.
The handset or clinician app is unresponsive	If the handset or the clinician app is unresponsive: 1. Turn the handset Off, then turn the power back On and make sure the handset is charged. 2. If that doesn't resolve the issue, use another handset, if available 3. If the issue persists, contact Medtronic support.
I can't find my downloaded report	Downloaded reports can be found in the Reports folder on the handset. You can access a downloaded report by navigating to the Reports folder on the handset or by connecting the handset to a PC and navigating to the Reports folder through the PC. See "Accessing a Downloaded Report" on page 57 for more information. Note: Generated reports will remain on the handset for 24 hours before being deleted.

Scenario	Explanations and possible solutions
I want to unpair or disconnect the neurostimulator from the apps.	Disconnecting or unpairing a connected neurostimulator is done using the Manage Devices feature. To disconnect: 1. Tap  to open a list of options, and select Manage Devices. Note: the Manage Devices feature can also be accessed by tapping  on the Menu screen. 2. Tap  if you see this icon. If you do not see the icon, proceed to step 3. 3. Tap X. 4. Tap REMOVE DEVICE to confirm that you want to unpair or disconnect the device from the clinician and patient apps. Note: Removing the device from the clinician app will disconnect the device from the patient therapy app. All data will also be erased from the app.

Table 8. Clinician app system notifications and solutions

App notification	Explanations and solutions
Permission required	Permission is needed to use the clinician app. 1. Tap EXIT. 2. Go to settings and select the clinician app. 3. Tap PERMISSIONS. 4. Slide from Off to On. 5. If the problem persists, contact Medtronic support.

System error	The system has encountered an unexpected problem. 1. Restart the clinician app. 2. If the problem persists, contact Medtronic support.
Low battery	The neurostimulator battery is low and therapy will be unavailable soon. Recharge the neurostimulator. (Model 97810 InterStim Micro neurostimulator) Replace the neurostimulator. (Models 3023 and 3058 neurostimulators)
Neurostimulator End of Service (Models 3023 and 3058 neurostimulators only)	The neurostimulator has reached end of service. The neurostimulator will no longer be able to provide therapy. You will need to replace the neurostimulator to continue therapy.
Neurostimulator Elective Replacement Indicator (ERI) (Model 97810 InterStim Micro rechargeable neurostimulator only)	The rechargeable neurostimulator battery is at or nearing the end of its service life. Consider replacing the neurostimulator.
Communicator connection lost	The communicator may be turned off, out of range, or the battery may be depleted. 1. Ensure the communicator is on. Move the communicator closer to the handset, and attempt to reconnect following the instructions presented on the screen. Note: If you are connecting to the communicator for the first time, see, "Using the Communicator to Pair the Handset and Neurostimulator" on page 16 for instructions.

	2. If the issue persists, the communicator or handset may need to be charged. Recharge the battery using the charger that came with the product. 3. If you are still experiencing issues, contact Medtronic support.
Communicator low battery	The communicator battery is below 25% and needs to be charged. The communicator battery LED indicator will be yellow when the battery is low. Charge the communicator using the charger that came in the product package. Note: If the communicator battery is below 10%, it will need to be charged above 10% before use.
Invalid Cycle time - Cycling error (Models 3023 and 3058 neurostimulators only)	The clinician app detects you are attempting to set the Cycling time to be shorter than the SoftStart/Stop feature time. Ensure the Cycling time is set to longer than the SoftStart/Stop feature time or tap CHANGE SOFTSTART/STOP FEATURE TIME to be shorter than the Cycling time. The clinician app detects you are attempting to set the Cycling time to one (1) second or less when the SoftStart/Stop feature is enabled. Ensure the Cycling time is set to greater than one (1) second when SoftStart/Stop feature is enabled or tap DISABLE SOFTSTART/STOP FEATURE to disable the SoftStart/Stop feature.
Invalid SoftStart/Stop Time-SoftStart/Stop Feature Error (Models 3023 and 3058 neurostimulators only)	You are attempting to set the SoftStart/Stop time longer than the Cycling time when Cycling is enabled. Ensure the SoftStart/Stop time is shorter than the Cycling time or tap CHANGE CYCLING TIME to change the Cycling On/Off time to the same value as SoftStart/Stop time.

	Ensure the future settings follow these rules: 1. Cycling On Time must be greater than or equal to SoftStart/Stop Time. 2. Cycling Off Time must be greater than or equal to SoftStart/Stop Time.
Data lost (Model 97810 InterStim Micro rechargeable neurostimulator)	Data has been lost. 1. Check rechargeable neurostimulator configuration and therapy programs. 2. If the problem persists, contact Medtronic support.
Data lost (Models 3023 and 3058 neurostimulators)	The neurostimulator has undergone a Power On Reset (POR). No therapy is available until the information in the neurostimulator is restored. Restart the clinician app. You will need to reconfigure therapy settings for the neurostimulator.
Power on Reset detected (Model 97810 InterStim Micro rechargeable neurostimulator only)	Device reset has occurred. Please check the neurostimulator battery level. 1. Check the rechargeable neurostimulator battery level and recharge the battery if needed. 2. Go to therapy programming to turn therapy on if needed. 3. If the problem persists, contact Medtronic support.
Invalid key	Your password key has expired. Refer to "I have forgotten my password" on page 67 and follow the steps shown in that section.
An invalid neurostimulator serial number is displayed on the Connection screen (e.g. NJY/NBV000000) (Models 3023 and 3058 neurostimulators only)	You are attempting to connect a neurostimulator that has experienced a Power On Reset (POR). 1. If the correct serial number is displayed, select the correct serial number, and tap CONTINUE. If the correct serial number is not displayed, skip to step 3.

	2. Select the correct serial number on the Connection screen, and tap NEXT to connect to the neurostimulator. 3. If the correct serial number is not displayed on the Connection screen, select the invalid neurostimulator serial number (e.g. NJY/NBV000000), and tap CONTINUE. 4. Enter the serial number of the neurostimulator and tap NEXT to connect to the neurostimulator. Notes: Once the neurostimulator is paired, the entered serial number will be displayed on the Connection screen for all subsequent connections. If you opt to manually enter the serial number, all stored data on the handset will be deleted.
Update required (Model 97810 InterStim Micro rechargeable neurostimulator only)	Updating the device will rewrite all settings. Click UPDATE to update the rechargeable neurostimulator, or END SESSION to not update.
Therapy state changed (Models 3023 and 3058 neurostimulators only)	Therapy output was changed by another programming device. Check therapy output. Use one handset and communicator at a time to make adjustments to the neurostimulator.
Maximum settings reached (Model 97810 InterStim Micro rechargeable neurostimulator only)	The neurostimulator battery is unable to produce the levels of energy required for the current stimulation settings. See "Maximum Settings Reached (Model 97810 only)" on page 73.

Maximum Settings Reached (Model 97810 only)

Considerations for troubleshooting

If a **Maximum Settings Reached** message appears during a programming session, the following are considerations for troubleshooting to return the system back to its expected functionality:

1. Take an impedance measurement to see if there is an open or a short circuit in the system.
2. If the **Maximum Settings Reached** message appears after a parameter has been changed, consider lowering that parameter.
3. Consider changing the electrodes that are currently providing stimulation.
4. If the battery level is low, consider recharging the neurostimulator.
5. If you are still experiencing issues, contact Medtronic support.

Clinician Reset and Recharge (Model 97810 only)

Perform Clinician Reset and Recharge if the clinician app or the recharger is unable to communicate with the neurostimulator, or if troubleshooting indicates a clinician reset is needed. Clinician Reset and Recharge will use the recharger to reset the neurostimulator.

Note: If the volume on the recharger is turned off, you will be unable to hear the audio feedback in the following steps. See the manual that came with the recharger for instructions to change the recharger volume.

Performing Clinician Reset and Recharge

1. Connect the charging dock to a power supply. Place the

recharger on the dock as if to charge the recharger.

2. Palpate the patient's skin to determine the approximate location of the neurostimulator.
3. With the recharger on the dock, enter Clinician Reset and Recharge by using the button press sequence: short, short, long. You will see the battery indicator fill up and hear a beep with each button press.

Note: A short button press is 1-2 seconds. A long button press is 4-6 seconds.

4. The light on the power button indicator will show a spinning amber light and the recharger will beep.
5. Pick up the recharger from the dock and hold it over the patient's implanted neurostimulator. Within 2 minutes, the power button indicator will show a flashing amber light and beeping will get faster. This means the recharger is preparing to reset the neurostimulator.
6. Keep the recharger over the implant after the beeping stops. The reset will start after the beeping ends and will take up to 1 minute to complete. The recharger indicator light will change to a slow pulse when the reset has completed.
7. After the reset, the recharger transitions to a special recharging session that may take up to 30 minutes if the neurostimulator battery needs to be recovered. Keep the recharger over the implant until it enters a normal recharging session.

Note: The recharger indicator light will turn green when the implanted neurostimulator enters a normal recharging session.

8. Check the therapy status after performing Clinician Reset and Recharge and turn therapy on if necessary.

For instructions on performing a normal recharging session on the Model 97810 InterStim Micro neurostimulator battery, see the manual provided with the recharger.

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
- Patient therapy app for sacral neuromodulation therapy - Patient User Guide for InterStim Systems

Ashley Lylloff
Pr. Regulatory Affairs Specialist
Medtronic Inc.

Date: 09 January 2023

Notary

Kimberly J. Horsch
January 10, 2023

A52200/
A520

PATIENT THERAPY APP
FOR SACRAL NEUROMODULATION
THERAPY

Patient User Guide for InterStim™
Systems

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The BLUETOOTH® word and logos are registered trademarks owned by Bluetooth SIG, Inc., and any use of such marks by Medtronic, Inc., is under license.

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Explanation of symbols

Manufacturer

Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts.

Authorized representative in the European community

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Importer

TABLE OF CONTENTS

HOW TO USE THIS GUIDE	7
Contact Information	8
How to Contact Medtronic	8
Important Tips.....	9
 GLOSSARY	 11
 GENERAL WARNINGS & PRECAUTIONS	 14
Warnings	14
Precautions	14
 ABOUT THE PATIENT THERAPY APP	 15
Intended use	15
How Does the App Work?	15
How do I use the patient therapy app on the handset?	16
Accessing the Patient Therapy App	17
Patient Therapy App Tutorial.....	18
Home Screen.....	18
Accessing system information	19
 USING THE PATIENT THERAPY APP	 21
Main Features	21
Using the Communicator to Connect your Handset and Neurostimulator	21
Adjusting Therapy Amplitude	24
Turning Therapy On or Off.....	24

Switching Programs.....	25
Checking Battery Level.....	25
InterStim Micro rechargeable neurostimulator battery	26
Table 1. Battery level descriptions (InterStim Micro rechargeable neurostimulator only)	27
InterStim and InterStim II battery level.....	27
Table 2. Battery level descriptions (InterStim and InterStim II only)	27
 MRI EXAMINATIONS	28
Preparing for Your MRI Appointment	28
Activating MRI Mode Using Your Handset ...	30
Table 3. MRI mode eligibility displays	32
Turning Therapy On After Your Scan.....	33
Visiting Your Clinician to Activate MRI Mode	33
 ADDITIONAL INFORMATION	35
Data Security and Network Connectivity	35
Software Updates	36
Safety and Technical Checks	36
Cleaning and Care	36
Disposal.....	37
 TROUBLESHOOTING	38
Therapy (stimulation) Problems.....	38
Table 4. Possible scenarios and solutions for therapy problems	38
App Notifications	47
Table 5. Important app notifications.....	48

Other Troubleshooting Scenarios..... 52

Table 6. Other troubleshooting scenarios and solutions53

HOW TO USE THIS GUIDE

Welcome to your guide to the patient therapy app for Medtronic Bladder or Bowel Control Therapy delivered by the InterStim™ system (also called sacral neuromodulation therapy). We created this guide to help you understand how the patient therapy application (app) on your handset works. Please read the entire guide before using the app to adjust your therapy. Ask your clinician to explain anything that is unclear.

- Information about how to contact Medtronic is on page 8 and the back cover.
- A glossary is provided in this guide to describe terms relating to your therapy that may be unfamiliar to you. Terms that are bold throughout this guide are defined in the glossary. Glossary terms and definitions begin on page 11.
- Important safety information regarding your therapy begins on page 14.
- The "USING THE PATIENT THERAPY APP" section provides instructions on how to use the app to adjust your therapy after your InterStim device has been implanted. This section of the guide includes information on turning your therapy on and off, adjusting your therapy settings, checking battery level, and changing programs (if applicable). See page 21.
- The "MRI EXAMINATIONS" section includes instructions for using the app to prepare for an MRI scan. See page 28.
- The "TROUBLESHOOTING" section provides explanations and solutions for

various troubleshooting scenarios you may encounter while using your handset. This includes explanations on important patient therapy app notifications, and who to contact if your handset is lost or broken. "TROUBLESHOOTING" starts on page 38.

Contact Information

Be sure to write down important contact information in the following fields, so it's easy to find in the future:

Clinician Name	Clinician Telephone Number
Hospital Name	Hospital Telephone Number

How to Contact Medtronic

Medtronic is available to answer any technical or troubleshooting device related questions or concerns you may have about your InterStim system components. Have the name and telephone number of your **clinician** at hand if you have any questions or problems. See the Medtronic contacts on the back cover of this manual.

! USA For assistance in the US, call 1-800-510-6735. Support is available Monday through Friday from 8:00AM to 5:00PM (central time).

If you lose your identification card, contact:

USA

Medtronic Inc., Patient Registration Services
Mail Stop SLK35
800 53rd Avenue NE
Minneapolis, MN 55421-1200
(1-800-551-5544)

Europe, Africa, Middle East, and Asia-Pacific Countries

Refer to the Medtronic contacts on the back cover of this manual.

Important Tips

The following includes some things to keep in mind while using sacral neuromodulation therapy:

- Always carry your handset and **communicator**. The patient therapy app on the handset is the only way to adjust or turn therapy on and off. Also, it is very important that you bring your handset and communicator to all follow-up appointments, including appointments with other clinicians. During certain procedures, you may need to turn your therapy off.
- Bring this guide to all your appointments and share this information with clinicians, dentists, MRI operators, and other professionals who may not be familiar with the InterStim system patient therapy app. When necessary, you should also bring your Patient Therapy Guide for InterStim systems; it contains important information

about sacral neuromodulation therapy that your clinicians should be aware of.

- Refer to "GENERAL WARNINGS & PRECAUTIONS" on page 14 and throughout the manual for information on what you should and should not do while using the InterStim system.

Notes:

For more information about the communicator, see the TM90 Communicator Instructions for Use.

For more information about the recharger (for use with the Model 97810 InterStim Micro rechargeable neurostimulator only) see the WR9220 Recharger Patient User Guide.

GLOSSARY

GLOSSARY

The following is a list of glossary terms that may be unfamiliar to you. The words that appear in **bold** in this guide are defined here for your convenience.

Caution	A statement that describes an action or situation which could harm the patient or damage the device.
Clinician	The term clinician is used in this guide to refer to a specialized medical professional such as a doctor, nurse, medical technician, or specialist.
Communicator	A communication device that includes a radio-frequency technology that allows the neurostimulator and the handset to communicate with each other.
Electromagnetic interference	A strong field of energy near electrical or magnetic devices that could prevent the neurostimulator from functioning properly, damage the neurostimulator, or even injure you.
Handset	A hand-held electronic device used by the patient to control therapy settings and to turn the neurostimulator on and off.

GLOSSARY

InterStim system	An implanted system that sends mild electrical pulses to sacral nerves in order to treat your symptoms.
Lead	A thin wire with protective coating that has metal electrodes on one end and a connector on the other.
Patient therapy app	An app used to modify sacral neuromodulation therapy.
Neurostimulator	The power source of a neurostimulation system. It contains the battery and electronics that control the stimulation.
Pair (concerning a device or devices)	The process of wirelessly connecting in order to send and receive information between the patient therapy app, the communicator, and the implanted neurostimulator.
Precaution	See Caution.
Program	A specific combination of stimulation settings assigned to deliver therapy to a specific site.
Sacral nerves	Nerves located near the base of the tailbone just above the cleft of the buttocks. Sacral nerves control the bladder, bowel, and pelvic organs.

11/25/2017

Stimulation	The delivery of electrical pulses to a specific site.
Therapy	Treatment of a disease or condition. When neuromodulation therapy is prescribed, a neurostimulation system is used to deliver therapy (stimulation) to a specific site.
Therapy amplitude	The strength or intensity of stimulation measured in Volts or Milliamps.
Therapy settings	Refers to the features that define the therapy. The clinician programs all therapy settings. You can adjust some of your therapy settings within clinician-defined limits.
Warning	A statement that describes an action or situation which could seriously harm the patient.

GENERAL WARNINGS & PRECAUTIONS

GENERAL WARNINGS & PRECAUTIONS

Warnings

 **WARNING: Wound contact** — Do not use the communicator on an unhealed wound. The communicator is not sterile and contact with the wound may cause an infection. Keep a sterile bandage or barrier between the wound and the device.

Precautions

 **CAUTION: Communication interference from EMI sources** – When using your handset and communicator to communicate with your **neurostimulator**, move away from equipment that may generate **electromagnetic interference (EMI)** or turn off the likely source of EMI. EMI may disrupt communication between the handset, communicator, and neurostimulator. Examples of EMI sources are computer monitors, and motorized wheelchairs. For more information about EMI, refer to the Patient Therapy Guide for InterStim systems.

 **CAUTION:** Contact your clinician if you experience unintended changes in therapy.

Other **warnings** and **precautions** are found throughout the manual.

For warnings, precautions, and adverse events related to sacral neuromodulation therapy, refer to the Patient Therapy Guide for InterStim systems. If you did not receive the Patient Therapy Guide, contact your clinician.

ABOUT THE PATIENT THERAPY APP

ABOUT THE PATIENT THERAPY APP

The **patient therapy app** is a software application (app) for Medtronic Bladder or Bowel Control Therapy delivered by the InterStim system. The app is intended to be used for controlling and monitoring your therapy. It is installed on and used with the **handset** to adjust your therapy.

Throughout this manual, "patient therapy app" refers to both the Micro My Therapy app used with the InterStim Micro SureScan™ MRI rechargeable neurostimulator (Model 97810), and the My Therapy app used with InterStim (Model 3023) and InterStim II (Model 3058) neurostimulators.

Intended use

The patient therapy application provides the user interface for patients to adjust therapy and obtain system information from an implantable neurostimulator as part of a neurostimulation system for sacral neuromodulation therapy.

How Does the App Work?

Your clinician used the clinician app on your handset to program your **therapy settings**. Your therapy is provided through the implanted neurostimulator (INS). Your clinician uses programs to control the delivery of your therapy. If needed, your clinician will use the clinician app to change your therapy settings. You will use the patient therapy app to:

- Turn your therapy on and off, page 24

ABOUT THE PATIENT THERAPY APP

- Adjust therapy amplitude within the limits set by your clinician, page 24
- Change **programs** (if available), page 25
- Activate MRI mode, page 30

The patient therapy app on the handset uses a communicator to communicate wirelessly with the implanted neurostimulator (INS) by sending signals to and receiving signals from the neurostimulator. To send and receive the signals, the communicator must be connected or **paired** to the handset as well as the neurostimulator. Your neurostimulator only accepts programming from the clinician app or the patient therapy app. Other devices are not able to program your neurostimulator.

Your clinician will initially set up the handset according to your specific needs, after the implant procedure. However, if for any reason the communicator and/or handset has been disconnected or unpaired from the neurostimulator, or you received a replacement device(s), you will need to reconnect. For more information about connecting the communicator to the handset and implanted device, see "TROUBLESHOOTING" on page 38.

Note: Make sure to keep your handset and communicator charged and with you at all times, in the event that you need to adjust or turn your therapy off.

How do I use the patient therapy app on the handset?

Use the patient therapy app on the handset to adjust your therapy as directed by your clinician.

ABOUT THE PATIENT THERAPY APP

The handset has a touchscreen display that allows you to tap or select icons that are displayed on the screen. Use one finger when tapping icons. You can use the touchscreen in the following ways depending on the function or task:

- Short, light taps for icon selection
- Swipe left/right to turn your therapy on or off.

Note: None of the physical buttons on the handset itself will alter your therapy.

Instructions for using the app to adjust your therapy settings begin on page 24.

Accessing the Patient Therapy App

You can access your patient therapy app (Figure 1) by locating the appropriate icon on the main screen of the handset. Tap the icon to open the app and begin using it.

The first time you or your clinician opens the app, it will request permission for data storage. Tap **OK** to continue. Then tap **ALLOW** for the next messages.

Figure 1. Micro My Therapy app icon (left) for the InterStim Micro rechargeable neurostimulator and My Therapy app icon (right) for InterStim and InterStim II

ABOUT THE PATIENT THERAPY APP

Note: Ask your clinician if you are unsure whether your neurostimulator uses the Micro My Therapy app or the My Therapy app.

Patient Therapy App Tutorial

The app includes a brief tutorial that will help you understand its basic functions. When you open the app for the first time you will be presented with this tutorial. If you would like to view the tutorial again in the future, you can do so at any time.

To view the tutorial:

Tap  to open the menu, and select **Tutorial** from the list of options.

Home Screen

The Home screen is the screen where you will adjust your therapy settings. Figure 2 provides an overview of the Home screen.

Note: The Home screen may look slightly different from the image shown depending on how your clinician set up the handset.

ABOUT THE PATIENT THERAPY APP

Figure 2. Home screen

Accessing system information

Use the patient therapy app to access neurostimulator system information as well as check the battery level of the neurostimulator, and prepare for MRI. (Figure 3).

ABOUT THE PATIENT THERAPY APP

Figure 3. Accessing system information

To view system information:

Tap  in the corner of the screen and select from the following:

- **About** to see general information about your system
- **Battery** for battery level (See “Checking Battery Level” on page 25 for details.)
- **MRI** to prepare for MRI (See “MRI EXAMINATIONS” on page 28 for important information.)
- **Tutorial** to review basic functions

USING THE PATIENT THERAPY APP

USING THE PATIENT THERAPY APP

Your clinician will provide complete guidelines about when you may want to adjust your therapy. The following sections provide instructions for using the app to modify your therapy.

Main Features

The patient therapy app has a few basic features that allow you to manage certain aspects of your therapy. Using the app, you can perform the following tasks:

- Adjust **therapy amplitude** (stimulation), page 24
- Turn your therapy On or Off, page 24
- Switch between **programs** (if available), page 25
- Access the app tutorial, page 18
- Access system status and neurostimulator battery level information, page 25
- Activate MRI mode, page 30

Using the Communicator to Connect your Handset and Neurostimulator

The patient therapy app will need to communicate with your implanted neurostimulator. Use the following procedure to use the communicator to connect your handset and implanted neurostimulator before adjusting your therapy.

USING THE PATIENT THERAPY APP

To connect your handset and neurostimulator using the communicator:

1. Press  on the communicator to turn the communicator on. The communicator will attempt to connect to the handset (Figure 4). The blue LED indicator on the communicator will continuously blink to indicate it is on and in discovery mode.

Note: For details about button and functions of the communicator, see the TM90 Communicator Instructions for Use.


Figure 4. Connection screen

2. Open the patient therapy app on the handset.
3. Once the communicator has successfully connected to the handset, the blue LED indicator on the communicator will be

USING THE PATIENT THERAPY APP

solid and no longer blinking. Place the communicator over the neurostimulator site (Figure 5), and tap **FIND DEVICE** on the handset.

Note: If the communicator fails to connect to the neurostimulator, readjust the location of the communicator over the neurostimulator and tap **RETRY**.

Figure 5. Connecting the communicator to the implanted neurostimulator

4. Once the communicator has successfully connected to the implanted neurostimulator, you will see the tutorial (if enabled) or the Home screen where you can adjust your therapy as directed by your clinician.

USING THE PATIENT THERAPY APP

Adjusting Therapy Amplitude

The patient therapy app allows you to adjust your therapy amplitude to your comfort level as your needs change.

To adjust your therapy amplitude:

1. Place the communicator over the neurostimulator site.
2. On the Home screen, tap the **arrows** to increase or decrease therapy amplitude as directed by your clinician.

Turning Therapy On or Off

In the patient therapy app you have the ability to turn therapy on or off.

To turn therapy on:

1. Place the communicator over the neurostimulator site.
2. On the Home screen, swipe the **On/Off** switch from Off to **On**.

Note: You can also turn therapy on by tapping **^** on the Home screen.

To turn therapy off:

1. Place the communicator over the neurostimulator site.
2. On the Home screen, swipe the **On/Off** switch from On to **Off**.
3. A pop-up will appear confirming that therapy has been turned off. Tap **OK** to acknowledge therapy has been turned off.

USING THE PATIENT THERAPY APP

Note: You can also turn therapy off by tapping  on the Home screen until the therapy amplitude is set to -- (off).

Switching Programs

Programs are created by your clinician to provide you with the best therapy results. Each program has different therapy settings. If multiple programs are made available to you, you have the ability to switch between the available programs as directed by your clinician.

Note: If your clinician set up your handset so that only one program is available to you, switching programs is not an option.

To switch between programs:

1. Place the communicator over the neurostimulator site.
2. Tap the **PROGRAM** button at the bottom of the Home screen.
3. Select the desired program from the list, and tap **DONE**.
4. Tap the **arrows** to adjust the therapy amplitude to the desired setting. See "Adjusting Therapy Amplitude" on page 24 for more information.

Checking Battery Level

To check neurostimulator battery level:

Tap  in the corner of the Home screen and select **Battery**.

USING THE PATIENT THERAPY APP

InterStim Micro rechargeable neurostimulator battery

⚠ CAUTION: Recharge a rechargeable neurostimulator as needed. A fully depleted neurostimulator battery will result in loss of therapy.

Be sure to check the battery level of the InterStim Micro rechargeable neurostimulator regularly. Figure 6 shows an example of the battery screen. Table 1 displays the battery level descriptions.


Figure 6. Battery screen (InterStim Micro rechargeable neurostimulator only)

USING THE PATIENT THERAPY APP

**Table 1. Battery level descriptions
(InterStim Micro rechargeable
neurostimulator only)**

Status	Battery level
Status: 30%-100% 	Battery Level: OK
Status: 20% 	Battery Level: LOW Recharge the neurostimulator.
Status: 10% 	Battery Level: CRITICAL Recharge the neurostimulator.

Note: For instructions on recharging the InterStim Micro neurostimulator battery, see the WR9220 Recharger Patient User Guide.

InterStim and InterStim II battery level

Table 2 displays the battery level descriptions for InterStim and InterStim II neurostimulators.

**Table 2. Battery level descriptions
(InterStim and InterStim II only)**

Status	Battery level
Status: 	Battery Level: OK
Status: 	Battery Level: LOW

CONFIDENTIAL

MRI EXAMINATIONS

Magnetic Resonance (MR) Conditional:

Depending on what kind of neurostimulation system components that you have implanted, you may be eligible for one of the following types of magnetic resonance imaging (MRI) scans:

- MRI scans of the head only (this is referred to as head-scan only eligible).
- MRI scans of any part of your body (this is referred to as full-body scan eligible).

Note: Always consult the clinician managing your InterStim system before having an MRI scan.

Magnetic Resonance (MR) Unsafe: External devices such as your handset or communicator are unsafe in the MR environment. Do not take any external devices into the MRI scanner (magnet) room.

Preparing for Your MRI Appointment

WARNING: Do not undergo an MRI scan without first determining if you are eligible and consulting with your clinician. Inform the MRI clinician that you have an implanted device, and activate MRI mode before your MRI scan, while you are still outside the scanner (magnet) room. Failure to do so may result in serious or permanent harm to you or damage to your implanted system.

CAUTION: If you have a rechargeable neurostimulator, ensure the battery level of the neurostimulator is at least 30% before undergoing an MRI scan. You are not able to receive an MRI

scan if your rechargeable neurostimulator is not sufficiently charged.

△ CAUTION: Deactivate MRI mode when the MRI scan is complete and you are outside the MRI scanner (magnet) room. Failure to deactivate MRI mode and return therapy amplitude to the previous level may result in a return of symptoms.

Before Your MRI Scan:

1. Consult the clinician managing your InterStim system before having an MRI scan.
2. If you have a rechargeable neurostimulator, ensure the battery level of the neurostimulator is at least 30% before undergoing an MRI scan.
3. Bring your handset, communicator, and Medtronic patient ID card to every MRI appointment, and inform the MRI clinician that you have an implanted device.
4. If the clinician managing your InterStim system gave you an MRI scan eligibility sheet, bring it to the MRI appointment along with your handset and communicator.
5. Activate MRI mode before your MRI scan, while you are still outside the scanner (magnet) room, using the following instructions. MRI mode turns therapy off and prepares your InterStim system for scanning.

Notes:

Inform the MRI clinician if you think you have a fever. Having a fever may affect your ability to have an MRI scan.

To obtain the MRI guidelines, the MRI clinician can refer to the contact information on the back

cover of this manual, or the MRI clinician can go to www.medtronic.com/mri.

Activating MRI Mode Using Your Handset

The handset must be used to place the InterStim system in MRI mode before an MRI scan. Ensure you are outside of the MRI scanner (magnet) room before proceeding with the following steps. When you use the handset to place your InterStim system in MRI mode, scan eligibility icons will appear on the screen (See Table 3 on page 32). Show that screen to the MRI clinician.

Note: During the MRI scan, leave the handset outside the scanner (magnet) room with the InterStim system in MRI mode. Do not tap **DEACTIVATE**. Therapy must remain off.

1. Press ☺ on the communicator to turn the communicator on. The communicator will attempt to connect to your handset, but cannot do so until the patient therapy app is launched. The blue LED light on the communicator will continuously blink to indicate it is on and in discovery mode.
2. Open the patient therapy app on the handset to initiate the connection process.
3. Once the communicator has successfully connected to the handset, the blue LED light on the communicator will be solid and no longer blinking. Place the communicator over the implanted neurostimulator and tap **FIND DEVICE** on the handset.

Note: If the communicator fails to connect to the neurostimulator, readjust the location of

the communicator over the neurostimulator,
and tap **RETRY** on the handset.

4. Once the communicator has successfully connected to the implanted neurostimulator, you will be taken to the Home screen.
5. Tap  in the corner of the Home screen and select **MRI**.

Note: If the handset screen says **Not Ready for MRI Scan** and **Not Eligible**, do not proceed with an MRI scan. Your system is not eligible. Contact your clinician.

6. Tap **ACTIVATE** to activate MRI mode.
7. One of the following messages will appear:
 - MR Conditional Full Body Scan Eligible
 - MR Conditional Head Scan Eligible with Transmit/Receive Head Coil
 - MRI eligibility cannot be determined

Note: If the handset screen says **MRI eligibility cannot be determined**, further assessment by the clinician who manages your InterStim system is necessary before scanning. Contact your clinician.

8. After activating MRI mode, show your handset to the MRI clinician with the screen displayed. Do not tap **DEACTIVATE** or turn therapy back **ON** until your MRI scan is complete and you are outside the scanner (magnet) room. See "Table 3. MRI mode eligibility displays" for descriptions.

Table 3. MRI mode eligibility displays

Screen image, text and description
 MR Conditional Full Body Scan Eligible
 MR Conditional Head Scan Eligible with Transmit/ Receive Head Coil
 MRI eligibility cannot be determined

Notes:

Your scan eligibility depends on several factors. Show your handset to the MRI clinician with the screen displayed. If the handset screen says **MRI eligibility cannot be determined** or any other unexpected message, do not proceed with an MRI scan. Contact the clinician who manages your InterStim system.

Do not take your handset, communicator, or any other patient control device into the MRI scanner (magnet) room.

Do not deactivate MRI mode until after the MRI scan is complete and you are outside the scanner (magnet) room.

Turning Therapy On After Your Scan

Deactivate MRI mode when the MRI scan is complete and you are outside the MRI scanner (magnet) room. You can deactivate MRI mode and turn therapy on using your handset, or you can return to the clinician managing your neurostimulation system to have therapy turned back on for you.

Complete the following steps to turn therapy on using your handset.

1. Connect to your neurostimulator as usual.
2. On the Home screen, tap  then select **MRI**.
3. Tap **DEACTIVATE** at the bottom of the screen.
4. Tap **DEACTIVATE** a second time when prompted.
5. Tap **YES** when prompted to return to previous therapy settings, if desired.

Note: Contact your clinician if you are unable to turn therapy on.

Visiting Your Clinician to Activate MRI Mode

Inform the clinician managing your InterStim system that an MRI examination has been prescribed for you and that you need MRI mode activated. Remember to bring your handset and communicator to every appointment, including the MRI appointment. Your clinician may place your system in MRI mode for you.

Keep your neurostimulation system in MRI mode during the MRI scan. The clinician may also give

you an MRI scan eligibility sheet to bring to your

MRI appointment. Give the MRI scan eligibility
sheet to the MRI clinician.

ADDITIONAL INFORMATION

The following sections provide additional information about your device, including information on data security and network connectivity, cleaning and care, and disposal information.

Data Security and Network Connectivity

The patient therapy app uses and stores data about your health and medical device. This information is protected by the app and handset. Personal files or data you choose to store or access elsewhere on the handset are not protected.

You are ultimately responsible for protecting the data that is stored on the handset. To help protect your data, Medtronic recommends the following:

- Secure your handset by disabling the Wi-Fi network connection when using the app.
- Use a trusted, secure Wi-Fi network when Internet access is needed.

Note: Connecting to a Wi-Fi network is not required to use your therapy, however network connectivity is needed for Medtronic app updates and communicator firmware updates.

- If your handset is lost or stolen, contact Medtronic right away using the contact information at the end of this guide.
- If you no longer need the handset to manage your therapy, see "Disposal" on page 37.

If you suspect the security of your handset has been compromised, stop using the app (if possible) and contact Medtronic to document and respond to the suspected incident.

Software Updates

You might be notified that you should update your software. Connect to Wi-Fi to receive the update. When the update is complete, you can turn Wi-Fi off again.

Safety and Technical Checks

Periodic safety and technical checks or periodic maintenance of the handset are not required. If the handset requires repair or is nonfunctional, contact Medtronic. Refer to the contact information at the end of this guide.

Cleaning and Care

Follow these guidelines to ensure that the handset and accessories function properly:

- Keep the handset and communicator out of the reach of children and pets.
- Use the handset and communicator only as explained to you by your clinician or as discussed in this manual and others provided to you by Medtronic.
- Follow all warnings and cautions in this guide, and the Patient Therapy Guide for InterStim systems provided by your clinician.
- Handle the handset and communicator with care. Do not drop, strike, or step on the device.

- Do not dismantle or tamper with the handset or communicator. Do not attempt to service them yourself.
- Clean the outside of the handset and communicator with a damp cloth when necessary.
- The handset and communicator are not waterproof. Do not allow moisture to get inside.

Disposal

Contact Medtronic when you no longer need your handset. Medtronic can help you delete your personal and health data, and dispose of the handset. Contact information for Medtronic can be found at the end of this manual.

TROUBLESHOOTING

The following sections contain information you may use to troubleshoot the patient therapy app, and answers to potential questions you might have while using it.

Therapy (stimulation) Problems

Table 4 provides troubleshooting scenarios you may encounter, as well as explanations and solutions to resolve problems if you are experiencing issues with your therapy.

Table 4. Possible scenarios and solutions for therapy problems

Scenario	Explanations and/or Solutions
Delayed therapy changes	You do not feel therapy right away after turning on the neurostimulator or you feel therapy after turning off the neurostimulator. Your clinician may have enabled SoftStart/Stop™ so that therapy starts and stops gradually. If SoftStart/Stop is enabled, it may take a few seconds for your neurostimulator to turn on and off. You may feel a residual effect after the neurostimulator is turned off.

Scenario	Explanations and/or Solutions
Intermittent therapy	You feel your therapy only some of the time. Your clinician may have programmed your neurostimulator to turn on and off at regular intervals. If therapy is uncomfortable or symptoms return, call your clinician.
Therapy area changes	You notice an unexpected change in where you feel the therapy (stimulation). Call your clinician.
Not receiving symptom relief	Therapy may be turned off, or the therapy (stimulation) might not be strong enough. 1. Check that your therapy is on, then increase the therapy amplitude as directed by your clinician. 2. If you have an InterStim Micro rechargeable neurostimulator, check battery level and recharge if necessary. 3. If these steps do not solve the problem, call your clinician.

Scenario	Explanations and/or Solutions
Uncomfortable or intolerable therapy (stimulation)	You are experiencing side effects from the therapy (stimulation). 1. Turn the therapy (stimulation) down or off. See "Adjusting Therapy Amplitude" on page 24 or "Turning Therapy On or Off" on page 24 for instructions. 2. If you are still experiencing uncomfortable therapy (stimulation), call your clinician.
Uncomfortable tapping sensation	You feel a tapping sensation that is too slow or too fast. 1. Turn the therapy (stimulation) down or off. See "Adjusting Therapy Amplitude" on page 24 or "Turning Therapy On or Off" on page 24 for instructions. 2. If you are still experiencing an uncomfortable tapping sensation, call your clinician.

TROUBLESHOOTING

Scenario	Explanations and/or Solutions
Therapy does not relieve your symptoms	1. Adjust therapy amplitude OR switch programs (if available). 2. If the issue persists, call your clinician.
You have unexpected changes in therapy	1. Turn your therapy off. 2. Turn your therapy on by pushing the up arrow once and slowly increasing therapy amplitude to the desired level OR switch programs (if available). 3. If the issue persists, call your clinician.
The handset, or the patient therapy app is unresponsive	Therapy may have stopped: 1. Turn off the power on the handset, then turn the power back on and reopen the patient therapy app. 2. If the issue persists, call your clinician.

www.medtronic.com

Scenario	Explanations and/or Solutions
You will be passing through a theft detector or security device You will be using potentially dangerous equipment You will be having a medical or dental procedure	Before engaging in any of these activities, see the Patient Therapy Guide for InterStim systems for details.
You received external defibrillation	Turn the therapy off and call your clinician.
The handset is unable to connect to the communicator	1. The communicator is not powered on. You will need to turn on the communicator to communicate with your implanted neurostimulator. Once the communicator is powered on, follow the connection steps (A) provided at the end of this section to reconnect your existing communicator.

Scenario	Explanations and/or Solutions
	2. The communicator is not charged or is in the process of charging. You will need to charge the communicator using the charger provided in your product package. Once charged, follow the connection steps (A) provided at the end of this section to reconnect your existing communicator. Note: The communicator cannot be used while charging. 3. The communicator is out of range. Make sure the communicator is near the handset. Once in range, follow the connection steps (A) provided at the end of this section to reconnect to your existing communicator.

Scenario	Explanations and/or Solutions
	4. You received a replacement communicator. You will need to connect the new communicator to the handset and the implanted device. Follow the connection steps (B) provided at the end of this section to connect to a new/replacement communicator device. 5. Your communicator needs to be updated. Follow the connection steps (C) provided at the end of this section to update your communicator. 6. If the issue persists, call your clinician.

Scenario	Explanations and/or Solutions
A. To reconnect with your neurostimulator using an existing communicator you have used before:	1. Make sure the communicator is turned on. Note: For details about button and functions of the communicator, see the TM90 Communicator Instructions For Use. 2. Open the patient therapy app on the handset. 3. When the communicator has successfully connected to the handset, place the communicator over the neurostimulator site, and tap FIND DEVICE on the handset. Once connected, you will be taken to the Home screen where you can adjust your therapy. 4. If the communicator you would like to connect to is not shown, tap RETRY. 5. If you are having trouble, call your clinician.

Scenario	Explanations and/or Solutions
B. To connect using a new or replacement communicator for the first time:	1. Make sure the communicator is turned on. Note: For details about button and functions of the communicator, see the TM90 Communicator Instructions For Use. 2. Open the patient therapy app on the handset. 3. After the app searches for and fails to find a communicator that has been paired previously, a screen will be displayed showing that the communicator is not found. Tap SWITCH COMMUNICATOR on the handset to connect. 4. Select the communicator serial number found on the back of the communicator that you would like to pair to and tap CONTINUE. The serial number of the communicator should match the serial number shown on the screen. Note: If the communicator you would like to connect to is not shown, tap RETRY. 5. If you are having trouble, call your clinician.

Scenario	Explanations and/or Solutions
C. To update your communicator: 1. Connect the handset to a Wi-Fi network. 2. Make sure the communicator is turned on. 3. Open the patient therapy app on the handset. 4. A screen will be displayed showing that an update is needed for the communicator. Tap UPDATE on the handset to update. Note: The update will take a few minutes to complete. 5. If you are having trouble, contact your clinician.	

App Notifications

Notifications can pop up on the handset while using the patient therapy app to notify you of problems or issues with the InterStim system. Table 5 provides important notifications along with explanations and solutions.

Table 5. Important app notifications

App notification	Explanations and/or Solutions
Communicator connection lost	1. Ensure the Bluetooth® on the handset is turned on. 2. Make sure the communicator is fully charged. 3. Reposition the communicator over the neurostimulator, and tap RETRY on the handset. 4. If the issue persists, call your clinician.
Device not responding	1. Reposition the communicator over the neurostimulator, and tap RETRY.
Internal device data lost or Data lost	Call your clinician.
Implanted neurostimulator (INS) low battery (InterStim Micro rechargeable neurostimulator only)	Your neurostimulator will need to be recharged soon. Recharge the neurostimulator.

App notification	Explanations and/or Solutions
Implanted neurostimulator (INS) low battery (InterStim and InterStim II neurostimulators only)	Your neurostimulator will need to be replaced soon (InterStim and InterStim II neurostimulators). Call your clinician.
Maximum Settings Reached (InterStim Micro rechargeable neurostimulator only) Note: In addition to the notification, an orange icon will appear.	Your neurostimulator has reached the limit of its therapy (stimulation) range. 1. Recharge the neurostimulator if needed. 2. Turn therapy amplitude down. 3. Swipe the On/Off switch to Off, then to On. 4. Slowly increase therapy amplitude. 5. If the notification is shown again, switch programs (if available) and adjust amplitude to the settings specified by your clinician. 6. If the issue persists, call your clinician.

Figure 10-1: Error Messages

App notification	Explanations and/or Solutions
Neurostimulator Elective Replacement Indicator (ERI) (InterStim Micro rechargeable neurostimulator only)	Your rechargeable neurostimulator is at or nearing the end of its service life. Call your clinician.
Neurostimulator End of Service (EOS) (InterStim and InterStim II neurostimulators only)	Your neurostimulator battery is depleted. Therapy can no longer be adjusted. Call your clinician to discuss replacement options.
System Error	The system has encountered an unexpected error. Call your clinician.
Device Communication Lost	Communication is lost or has been interrupted. 1. Close the app and try to reconnect. 2. If the issue persists, call your clinician.

App notification	Explanations and/or Solutions
Not Installed Properly	The patient therapy app was not installed properly Contact Medtronic support.
Communicator Low Battery	The communicator battery is low, and may be unable to communicate with the handset and/or the neurostimulator soon. Charge the communicator. Note: If the communicator battery level is below 10%, the battery will need to be charged to above 10% before use.
No Implant Date (InterStim Micro rechargeable neurostimulator only)	Implant date not recorded in handset. Call your clinician.
Update required (InterStim Micro rechargeable neurostimulator only)	Device must be updated to continue. Call your clinician.

Table 6: Troubleshooting Scenarios

App notification	Explanations and/or Solutions
Power on Reset detected (InterStim Micro rechargeable neurostimulator only)	POR has occurred. 1. Check neurostimulator battery level and recharge if needed. 2. From the Home screen, swipe the On/Off switch from Off to On. 3. If the issue persists, call your clinician.

Other Troubleshooting Scenarios

Table 6 provides additional troubleshooting scenarios and solutions.

Table 6. Other troubleshooting scenarios and solutions

Scenario	Explanations and/or Solutions
The handset has no power or has lost power	1. The handset battery is depleted. Recharge the handset using the charger. 2. The handset cannot be charged with the charger, or the handset may be damaged or malfunctioning. Contact Medtronic support.
The handset will not charge	1. The charger is disconnected from the handset. Connect the charger to the handset. 2. The charging cable or wall plug are not fully inserted. Ensure the charging cable and wall plug are fully inserted. 3. An incompatible charger is connected to the handset. Connect an appropriate charger to the handset. 4. The charger is defective. Call your clinician. 5. The handset is damaged or malfunctioning. Contact Medtronic support.

4.6" x 7.6" inches (116.84mm x 193.04mm)

Scenario	Explanations and/or Solutions
The handset is lost or broken	Contact Medtronic to order a replacement. Medtronic can also attempt to remotely delete your personal and health data from a missing handset. Contact information for Medtronic can be found on the back cover of this manual. You should also call your clinician to say that you have ordered a replacement.

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All Rights Reserved