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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SCENKER BIOLOGICAL TECHNOLOGY CO., LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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
Authorized
Signature:

ShaoYunjie

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(Date: Sep. 23, 2019)

证明 查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY

SCENKER Biological Technology
Co., Ltd.

CUI GE/Legal

representative

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September 12, 2019

Automated bacteriological system ЮЮЮЮЮ® LABSTAR 50
with the accessories

SERVICE MANUAL

Preface

We thank you for buying the Automated bacteriological system LABSTAR 50.

LABSTAR 50 provides the rapid growth of microorganisms and is designated for the clinical culture and detection of the pathogens in the blood and other aseptic body fluids. The system allows quickly and accurately detecting the pathogenic microorganisms in the blood of the patients with septicemia, bacteremia, etc. Besides, this device can also be used to detect the pathogenic microorganisms in the body cavities that are usually without bacteria (such as thoracic cavity, abdominal cavity, articular cavity, pericardial cavity and cerebrospinal cavity, etc.), that allows making the clinical diagnosis quickly and prescribing the effective treatment regimen.

This service manual contains the information necessary for the device operation and maintenance. The user should carefully follow the instructions described in the manual. It will allow extending the device service lifetime.

System delivery package includes the power-supply cord providing the safe operation and good ground of the device. Take care to ensure that the power-supply cord is directly connected to the triphase electrical outlet with the good ground. This ensures the device safe operation.

The System is recommended to be used in the laboratory that meets the requirements of Level-2 bio-safety laboratory. During the device operation, it is necessary to take the bio-safety controls. Used bottles must be sterilized by autoclaving or disinfected using other methods established by local regulations before disposal.

The System does not have the any parts that the user can replace on site. Do not disassemble the device body without any professional help.

Before starting the device operation, it is recommended to carefully read this manual. In case of any problems take place during the device operation or if you need the technical consultation, please address the manufacturer's authorized official in the Russian Federation directly.

The instrument will be presented under **ЮХОНА®** trademark on the Russian market.

The product information:

Production license No.: L.S.Y.J.X.S.C.X. No. 20110052

Manufacture date: Please, check out the label.

Product standard code: YZB/ LU 0278-2012 dated February 01, 2014.

The manual initial issue date: February 01, 2014

Revision date: December 01, 2014

The manual revision for the registration in Russian Federation:
v.2 dated September 12, 2019.

Scenker Biological Technology Co., Ltd.

Contents

1 Medical device data.....	5
1.1 Medical device name.....	5
1.2 Information about the Medical device manufacturer	5
1.3 Information about the Medical device manufacturer authorized official	5
2 Operation and description.....	5
2.1 Medical device overview.....	5
2.2 Technical details, parameters and characteristics of the medical device.....	6
2.3 The contents (contents of delivery) of the medical device.....	23
2.4 Medical device precaution measures and limits of use	23
2.5 The description of the package, marking and warning signs on the device	25
3 Medical device operation	28
3.1 Environmental requirements	28
3.2 Installation procedure	29
3.3 Setting of system parameters.....	33
3.4 Operation procedure	34
3.5 Culture media and application scope.....	46
3.6 How to collect the blood and other sterile biological fluid samples.....	46
3.7 Results explanation.....	49
4 Maintenance and repairs.....	50
4.1 The user-performed technical maintenance:.....	50
4.2 The common failures and treatment measures	51
5 Transportation and storage	52
5.1 Transportation	52
5.2 Storage.....	53
6 Warranty policy	53
7 The medical device utilization and elimination procedure.....	53
8 Appendix A	55
9 Appendix B.....	56
10 Appendix C.....	67
11 Appendix D	69
12 Appendix E.....	72

1 Medical device data

1.1 Medical device name

Automated bacteriological system **IOHOHa®** LABSTAR 50 with the accessories (hereinafter referred to as System).

1.2 Information about the Medical device manufacturer

SCENKER Biological Technology Co., Ltd., East End of Weier Road, Fenghuang Industrial Park, Liaocheng, Shandong Province, China, 252000, tel. (86) 800-860-8018

1.3 Information about the Medical device manufacturer authorized official

Concerning the quality and service of the Medical device “Automated bacteriological system **IOHOHa®** LABSTAR, design versions: **IOHOHa®** LABSTAR 50, **IOHOHa®** LABSTAR 100 with the accessories”, manufactured by SCENKER Biological Technology Co., Ltd., China, please contact the manufacturer’s authorized official in Russian Federation:

Triopromtrade, LLC,

The Russian Federation, 620028, the Sverdlovsk Region, Yekaterinburg, 28, Kirova St./43, Toledova St., Phone: +7 (343) 270-75-00, Fax: +7 (343) 246-91-56.
E-mail: info.triopromtrade@mail.ru

2 Operation and description

2.1 Medical device overview

2.1.1 Medical device designation specified by manufacturer

System is designated for the clinical culture and detection of the pathogens in the blood and other aseptic body fluids by the fluorescent method.

2.1.2 Scope of use and potential customers of the medical device

It is used in the medical and clinic diagnostic centers. The users are the qualified health care professionals engaged in the microbiological studies.

2.1.3 Indications for use

Only for in vitro diagnostics. The system allows quickly and accurately detecting the pathogenic microorganisms in the blood and other body fluids of the patients with septicemia and bacteremia, using the special culture bottles produced by SCENKER Biological Technology Co., Ltd.:

- Culture Medium Bottles **ЮHOHA®** for aerobe cultivation (color of label and removable cover is red, designation ■) – a standard culture medium bottles designed for adult patients not treated with antibiotics, the volume of the biological sample is from 5 ml to 10 ml;
- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for aerobe cultivation (color of label and a removable cover is blue, designation ■), culture medium bottles designed for adult patients treated with antibiotics, the volume of the biological sample is from 5 ml to 10 ml;
- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for anaerobe cultivation (color of label and removable cover is cadmium orange, designation ■), culture medium bottles designed for adult patients treated with antibiotics, the volume of the biological sample is from 2 ml to 10 ml;
- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for aerobe cultivation for children (the color of the label and removable cover is pink, designation ■), culture medium bottles designed for children, the volume of the biological sample is from 2 ml to 10 ml.

2.1.4 Contra indications:

Never use System with other manufacturers' culture bottles.

Never use System for the purposes not allowed by the manufacturer.

2.1.5 Precautions

System initial installation must be carried out by the specialist assigned by the manufacturer's authorized official in the Russian Federation for this purpose.

Before commencing the work with System, it is necessary to read the Service Manual to provide the equipment safe and efficient operation.

2.1.6 Risk class and applicable classification rules

The potential risk class of the medical device application according to the nomenclature classification of the medical devices is 2a (article 9.6 of the Order of the Minister of Healthcare of the Russian Federation as of 06/06/2012. No. 4H).

According to the Directive "On medical devices" issued by Chinese Medical Devices and Pharmaceuticals Control Department in 2002, the device is classified as "clinical facilities for testing and analysis 6840" and belongs to the second class of the medical devices.

2.2 Technical details, parameters and characteristics of the medical device

2.2.1 Description

This System is the automated system for continuous culture and analysis. The incubation of the culture bottles with their simultaneous shaking to stimulate the growth and monitoring the fluorescence change takes place in the main module of System without the personnel intervention.

The fluorescent method for microorganism detection is used in System. There is an optic sensitive element (polymer sensor) within the bottoms of the culture bottles which starts glowing as affected by the products of the bacteria metabolism, and also, the luminous intensity increases with the increase of the bacteria quantity. The system determines the microorganism growth depending on the nature of the luminance change.

The System allows quickly and accurately detecting the pathogenic microorganisms in the blood of the patients with septicemia, bacteremia, etc. Besides, this medical device can also be used to detect the pathogenic microorganisms in the body cavities that are usually without bacteria (such as thoracic cavity, abdominal cavity, articular cavity, pericardial cavity and cerebrospinal cavity, etc.), allowing making the clinical diagnosis quickly and prescribing the effective treatment regimen.

The Systems are designed for operation only with the special culture bottles containing the reagents necessary for the microorganism detection in the sample.

Each 10 minutes the system processes the signal received from the sensors. When detecting the positive samples, the visual notification system on the monitor is activated. More than 90% of the positive samples are detected within 24 hours.

In case no growth during the maximum fixed culture time, the system provides the negative result by the visual alert. After that the culture bottles can be removed from the device.

As the system operates with each culture bottle separately, the new samples determination can start in any moment.

There is a Barcode scanner for the convenience in operation.

The system has the following functions:

- The help function in case of the culture bottle incorrect installation.
- The automatic detection and calibration of the detected culture bottle position, and automatic quality control.
- System computer interface allows displaying the graphic information about the culture bottles position in the system under the condition of incubation.

- The provision is made for preparing the data reports on the testing at any time.

The system is adjustable and reliable in operation.

2.2.2 Medical device components, versions list


System **ЮHOHa®** LABSTAR 50 consists of (fig.1):

- Main module;
- Main module power cord;
- Computer with the preinstalled software;
- Data cable;
- Computer power cord.
- Barcode scanner (if necessary).

2.2.3 List of accessories

The accessories provided for System operation:

- Monitor with the power cord;
- Keyboard;
- Mouse pointing control;
- Signal wire;

*Fig. 1 – Automated bacteriological system **ЮHOHa®** LABSTAR 50 with the accessories
1 – Main module, 2 – Computer, 3 – Monitor, 4 – Keyboard, 5 – Mouse, 6 – Barcode scanner,
7 – Power indicator, 8 – Temperature indicator, 9 – Software*

2.2.4 General description of the main functional units, parts, body, functionality

2.2.4.1 Main module

System main module consists of:

- Revolute systems with the Positions for the bottles,
- Electric drive,
- Photoelectrical detection system,
- Temperature control system,
- Main control PWB.

LABSTAR 50 has the construction of the rotary turntable. The operation cycle (scan mode) of System occurs every 10 minutes. The analytical unit of the system computer creates the growth curve based on the measurement data, and its form allows making the conclusion of the presence or the absence of the pathogenic microorganisms in this sample.

Fig.2 - LABSTAR 50 System rotary turntable

The rotary turntable is the round flat surface with the diameter of 460 ± 1 mm and the thickness of 80 ± 1 mm. The flat surface has the holes – the positions to hold the bottles. The hole diameter is $36.5(+0.1-0.3)$ mm.

The holes for the bottles on LABSTAR 50 rotary turntable are divided into three round sections: inner (A), medium (B) and outer (C). The holding capacity is 50 bottles. The positions are numbered from 1 to 50: A01-A10, B01-B17, C01-C23. The positions are related to three sets of the highly sensitive fluorescent sensors.


Fig.3 –LABSTAR 50 main screen. Growth curve

The electric drive is the controlled EMS designed for the conversion of electrical energy into mechanical and vice versa, and the management of this process.


Fig. 4 – Gear-motor drive

The system consists of the gear-motor drive - the unit combining the functions of the reducing gear and the electric motor. The combined type five-phase step motor, with the AK-GB series reduction gear manufactured by Autonics Corporation, Korea, is used as the electric motor.

The drive is designed for the control of the step motor. The system uses the MD5-HF14 drive manufactured by Autonics Corporation, Korea.


Fig. 5 - The drive

In its turn, the drive receives the commands from the higher class controlling systems – PMC series controller manufactured by Autonics Corporation, Korea

The photoelectrical detection system is the detector block consisting of the sensor suite.


Fig. 6 - Photoelectrical detection system

The sensor consists of 4 generating tubes and 1 receiving photo element. Four symmetrically located red LEDs are used as the source. The photoelectric element is well centered towards the reference axes of four generating tubes. That provides more proportional light transmission. The light reflected from the bottle bottom, goes to the light cell receiver, where the optical signals are converted into the electrical signals. Conversely, these signals go through the units of amplifier and conversion of the analog signal into the discrete signal. The discrete signals are processed in the analytical unit of the system computer. The detector wave-length is: 630-650nm.


Fig.7 – The detection principle with the system usage

*1- transmitting tube, 2- polymeric sensor, 3-bottle with the medium,
4-reflected light, 5- receiving tube.*

The temperature control system is the system that provides the fixed temperature inside of System required for the correct culture. The system consists

of the heater, cooling fan, temperature detector and regulator device. The temperature detector transmits the value of the current temperature to the regulating device generating the control action based on the received information. Then the signal passes to the regulating device. Its task is to perform the functions to maintain the fixed temperature at the necessary level. In case the culture temperature falls outside the scope of the acceptable range, the signal is addressed to the controllers of the main control PWB controllers.

The main control PWB is the board that has the power check panel and the signal reception controllers.


Fig. 8 – Main control PWB


Fig. 9 - LABSTAR 50 System process functional diagram

2.2.4.2 Main module power cord

The power cord is designed for the main module power connection.

It has the following technical features:

Chinese Standard: CPCS-CCC

Model: PSB-10/ST3

Type I connector, 3 flat contacts.

Load-carrying capability: 10A/250B

Length, min.: 1 m


Fig. 10 - Main module power cord

2.2.4.3 The computer with the preinstalled software.

a. The computer has the following characteristics listed in Table 1

Table 1. Computer characteristics

Technical characteristics	LABSTAR 50
Operating system (OS)	Microsoft Windows Professional 10
OS license number	68536413
Processor type	Intel ^R Celeron ^R G1840
Number of cores	2
Maximum clock frequency	2.8 GHz
Memory cache	2 Mb
Socket	LGA 1150
Type	DDR3
Volume	2.00 Gb
Maximum RAM	4 Gb
Frequency	1,333 MHz
Hard disk (HDD)	500 Gb
Graphics card	Intel GMA HD 4000
Audio	Integrated
Power adapter capacity	AC220-240V ≥200W
Type/Quantity	USB 2.0 type A/4, USB 3.0 type All-WO/1, PS/2/2, RS232, RJ45, VGA
Disc drive type	DVD
Software	Preinstalled, LABSTAR 50, version 2.0,

	Copyright certificate of registration No. RZDZ 0251777
Overall dimension, mm	155x400x370

b. Software

The preinstalled PC software is the integral part of System and it works only together with all its parts. The software has the criteria that determine the signal processing when the strength of the fluorescence increases during the culture. The changes can be observed using the growth curve displayed on the monitor for each loaded position. The monitor also displays the information about the system status: temperature, time before the next testing, the total number of loaded bottles and culture data.

The physical and computer ability characteristics of the software are listed in Table 2

Software characteristics	Model LABSTAR 50
Name	LABSTAR 50
Version	V2.0
Platform	Intel ^R Celeron ^R G1840
Operating system	Microsoft Windows Professional 10
Software programming language	Delphi
Software data transmission communication interface	RS-232
Software data computing communication interface	Automatic
Software data storage communication interface	-
Software response at full load, s, max	1
Accident recovery time (complete software failure due to a technical or operational accident), h, max	2
Processing speed, Mbit /s	100
Recovery time, if virus programs are detected, h, max	2
Performance, Gb/s, min	1.5
Search time at full load, s, max	1
Interfaced systems	Ethernet(TCP/IP)
Protection means	Software account
Safety class	P "Minor injuries are possible"


Fig.11 - LABSTAR 50 screen with the system information

The monitor screen always displays the main interface:


Fig. 12 – LABSTAR 50 Main computer interface
1-Bottles position, 2- control panel, 3- growth curve,
4- system information

(1) – The position of the bottles on the turntable is always displayed, the empty positions have the white color, the loaded positions are yellow, the positive bottle positions are red and the negative bottle positions are green. The error positions are blue.

(2) - The control panel located in the screen footer consists of two buttons:

Open	For the main module door release
Install the bottle	For loading the bottles. The button is inactive when

	the door is closed.
Take out the bottle	For unloading the bottles. The button is inactive when the door is closed.
Data	For patient data input and request with the growth curve.
System setting	For obtaining the system information.
Close	For closing the program.

(3) – When passing the mouse to the loaded position, the diagram displays the growth curve. The data readout takes place each 10 minutes. The form of the curve reflects the presence or the absence of the microorganism growth in the loaded bottles.

By default, the range of Y-axis is 700-4,000, the range of X-axis is 0 - the maximum time of cultivating this bottles in terms of the days. If the bottle is defined as the positive one, the chart shows the “alarming point” - the detection point and the cultivating time before the bottle positive definition. The scale of X and Y axes can be changed using the right mouse button. The date and time of the bottle installation and the date and time of alerting are pointed out under the growth curve.

(4) – System information:


- Temperature – culture/incubation temperature;
- State of the system: “normal” at the normal temperature, “cooling” – if the temperature is over the limit, “heating” – if the temperature is below normal;
- Next testing – time before the next scanning/data readout;
- Operational state: “incubation” – culture process takes place, “scanning” – the data readout process takes place, “door open” – the main module is open for loading/unloading the bottles;
- The total number – the total number of the loaded bottles;
- In progress – the number of the bottles that are in the culture/incubation process;
- Positive – the number of the positive bottles;
- Negative – the number of the negative bottles.

More detailed information about the software functionality is provided in Article 3 “Medical device operation”.

2.2.4.4 Data cable

The data cable designed for connection of the main module to the computer.

Table 3. Technical characteristics of the data cable

Parameter	LABSTAR 50
Description	Multiway 1-3
Length, no less than, m	1.3
Load-carrying capability, A	0.5
Port	(RS232/37/37---37)
Photo	

2.2.4.5 Computer power cord

The computer power cord is designed for the computer power connection.

It has the following technical characteristics:

Chinese Standard: CPCS-CCC

Model: PSB-10/ST3

Type I connector, 3 flat contacts.

Load-carrying capability: 10A/250V

Length, min: 1 m


Fig. 13 - The computer power cord

2.2.4.6 Barcode scanner

The technical characteristics of the barcode scanner are listed in Table 4.

Table 4. Technical characteristics of the barcode scanner

Parameter	Characteristic
Model	LS-1881

Manufacturer	Guangzhou OWL Electronic Technology Co. , Ltd
Projector	Laser
Scanning frequency	150 scan/sec
Barcode contrast	≥30%
Wave-length	650 nm
Operational temperature	-45~45°C
Resolution	PCS 0.9 (0.1~0.625MM)
Focus	10mm~420mm
Indication (LED, sound etc.)	Double color LED
Protection level	IP54
Data transmission rate	200 characters per second
The mean of communication of the scanner with System;	RS232 serial data transmission
Readable codes	Single two-dimensional code

Fig. 14 - Barcode scanner

2.2.4.7 System parameters and characteristics

Table 5.

Index	LABSTAR 50
Bottles capacity	50
Culture temperature	35±1.5°C
Rotating speed	26 r/min
Operation cycle (scan mode)	every 10 minutes
Rotating mechanism type	Turntable
The system is designed for microorganisms culture and detection	+

The rotating and shaking of the bottles is provided	+
Automatic correction and quality control	+
The number of positions for the bottles	50
The availability of the graphic interface	+
The availability of the Barcode scanner	+
Fluorescent method (instrumental estimation), bottle bottom color change (manual estimation)	+
Incubation conditions: fixed temperature, rotating and vibration	+
A/C power source, V	220 ± 22
Current frequency, Hz	50 ± 1
Input power, max VA	450
Main module external dimensions, mm	$(640 \times 660 \times 650) \pm 10$
System main module weight, kg	98 ± 2
Efficiency*	10 samples a day
Time of cultivation	120 hours
Volume of biological sample	5-10 ml aerobic bottles and 2-10 ml of anaerobic and children

* The calculation is based on the full load of the analyzer and the time of cultivation

2.2.4.8 Basic requirements

a. Operating conditions

- Environment temperature: 10-30°C;
- General humidity: 45-85%;
- Atmosphere pressure: 86.0-106.0 kPa;
- Electric power conditions: 220V $\pm$ 22V, 50Hz $\pm$ 1Hz;
- Input power: $\leq$ 450VA;
- Away from the sources of the high electromagnetic disturbance and bright directed light;
- System should be dust-proof; the air should not contain the acids, salts and other corrosive gases;
- The proper grounding should be provided;
- It is necessary to control the vibration at the workplace.

b. Appearance

- The equipment should be flat, clean and of the uniform color, it should be without the scratches, blinks and mechanic damages.

- The equipment marking should be easy-to-read.

c. Controllability

The start button and the interconnecting wire should have the characteristics that facilitate comfortable usage, operation speed and reliability. The displayed information should be clear and constant.

d. Culture temperature

The temperature measurement accuracy and the oscillations should meet the following requirements:

- The culture temperature should remain within the range of $35 \pm 1.5^{\circ}\text{C}$.

- The temperature variations should not exceed 3.0°C .

e. System functions

- Continuous culture.

- Automated monitoring and determination of the positive and negative results.

- Immediate notification about the positive bottles detection.

- Alarm signaling in case the temperature exceeds the limits of the acceptable range.

f. Positive results repeatability

The results of the testing using the standard strains of the American typical culture collection (ATCC) should be positive (see Annex A).

g. Display functions

The condition of the bottle contents is determined as positive or negative. All the indicators are displayed on the system screen, regardless of the fact if the bottle with a culture corresponds to a detector particular position, and if the positive or negative results are obtained.

h. Barcode input and position

The system can recognize and place the bottle on the definite detector level according to the barcode.

i. Safety requirements

See the safety requirements in Appendix B,C, D.

2.2.5 The accessories overview, the contents and the functional capabilities

2.2.5.1 The monitor with the power cable

The main technical characteristics of the monitors are listed in Table 6.

Table 6. The monitor technical characteristics


Characteristics	LABSTAR 50
Brand, model	AOC, model E2070Swn
Photo	
Screen size	19.5'
Matrix type	TN
Screen format	16:9
Resolution	1,600x900
Brightness	200 kcd/m2
Contrast ratio	600:1
Response time	5 ms
Connector type	VGA
Power cable	Chinese Standard: CPCS-CCC Model: PSB-10/ST3 Type I connector, 3 flat contacts. Load-carrying capability: 10A/250B Length, min: 1.5 m

2.2.5.2 Keyboard

The computer keyboard is the device to input the data into the computer (input device). It is made as the set of the keys (buttons) arranged in the certain order.

Table 7. Keyboard technical characteristics

Characteristics	LABSTAR 50
Description	Standard classical desktop corded keyboard with 104 keys
Cable length, min, m	1.5
Connection type	PS/2

Photo	
-------	--

2.2.5.3 Mouse pointing control

Computer mouse is the pointing device for the cursor control and issuing the commands to the computer. The cursor control is made by moving the mouse over the table surface or mouse pad. The mouse buttons and wheel cause some actions, e. g.: selected object activation, webpage and e-document vertical scrolling.

Table 8. Computer mouse technical characteristics

Characteristics	LABSTAR 50
Description	Optical corded mouse
Cable length, min, m	1.5
Connection type	PS/2
Photo	

2.2.5.4 Signal wire


The low current signal wire is designed for connecting the computer with the monitor.

It has the following technical characteristics:

Length of min 1.3 m

Load-carrying capability 0.5 A

Connector: VGA (DE15F) — 15-contact subminiature DIX-connector connector for connecting the analogue monitors according to VGA graphic interface standard.

2.2.6 Sterility information

System is not sterile. Only the consumable bottles with the medium used with System, sterilized.

2.3 The contents (contents of delivery) of the medical device

- Main module – 1 pcs.
- Data cable - 1 pcs.
- Main module power cord – 1 pcs.
- The PC with the preinstalled software – 1 pcs.
- Computer power cord – 1 pcs.
- Barcode scanner – 1 pcs. (if necessary)
- Service manual – 1 pcs.
- Certificate – 1 pcs.

Accessories:

- Monitor with the power cable - 1 pcs.
- Cable - 1 pcs.
- Mouse pointing control – 1 pcs.
- Signal wire -1 pcs.

2.4 Medical device precaution measures and limits of use

2.4.1 Critical information

 It is not allowed installing the programs or files on the device computer, that are not related to its operation. In case of the computer failure or data loss for that reason, the manufacturer can carry out the maintenance only on the paying basis.

 During the operation of System the operator can be subjected to dangerous microorganisms, so it is recommended to install the device in laboratory premises that meet the requirements of biosafety level 2. The provisions should be made for the biological protection measures and it is necessary to fully comply with the requirements for the organization of the microbiological research laboratories. The used bottles and other medical wastes must be disposed of by

autoclaving sterilization or other methods determined by the local rules and regulations (see 7 The medical device utilization and elimination procedure).

There is the danger of infecting during testing the samples with System. As a result, the marking in the form of the Biohazard Mark is provided on the device.


It is prohibited to touch the electric power supply of the operating System.

The corresponding parts of the device have the “electric shock hazard” mark.

The device back panel has the mark “Caution! Please, refer to the service manual”.

In case of the noncompliance with the operating procedure established in this document, the normal operation of System can be disrupted. In this case the findings can be distorted, and the destruction of the equipment as well. The user and (or users) of the device hold responsibility of the consequences of this noncompliance.

System main module has the drive unit connected to the A/C power source with the voltage of 220 V and it can be a source of the external electromagnetic disturbance. (see Appendix D EMS)

2.4.2 FAQs

2.4.2.1 The pathogenic bacteria of the most frequent occurrence (>90%):

- Staphylococcus aureus
- Escherichia coli
- Enterobacteriaceae: Salmonella, Klebsiella pneumoniae etc.
- Pseudomonas aeruginosa
- Streptococcus pneumoniae
- Candida

2.4.2.2 The most commonly occurred infesting bacteria (3%-5%)

- Staphylococcus epidermidis
- Propionibacterium acnes
- Fusobacterium
- Diphtheria group
- Gram positive bacillus
- Acinetobacter

2.4.2.3 The reasons of microorganism low rate of growth, corrective measures

- Insufficient sample volume: within the scope of the sample volume (Adults: 5-10 ml, children: 2-10 ml) proportionally to the positive growth of the cultivated microorganisms.
- The sample collection within 24 hours is not enough: the bacteremia is remitted in 85% cases, several samples collection can assure the positive growth recognition.
- The blood sampling should be done within one locus every time.
- The inappropriate time of the blood sampling.
- Carrying out the antimicrobial therapy before or during the blood sampling.
- Storing the samples at the room temperature: it is not allowed storing the samples in the freezer, because the bacteria with the special nutritive demands die at the low temperatures.

2.5 The description of the package, marking and warning signs on the device

2.5.1 Package

System consumer package is the transport package at the same time.

The transport package should be made of woods (packing-box).

The packing-box overall dimensions:

Parameter	LABSTAR 50
External dimensions, mm	(760x790x790)±20
Weight, kg	24±2


Fig. 16 - LABSTAR 50 transport package

System is placed into the waterproof poly film and packed in the packing box that has to be furnished with the soft inserts between the device and package walls for the impact absorption.

The operational documentation is put into the separate poly bag and packed into the box.

The computer with the preinstalled software is packed and transported in the manufacturer's cardboard pack. The accessories are packed there as well: the keyboard and the mouse.

The monitor is transported in the manufacturer's cardboard pack.

2.5.2 The product marking

The LABSTAR 50 System product marking contains:

On the front panel of the main module (door):

- Manufacturer's logo;
- Name LABSTAR 50;
- Indicator and "power supply" caption;
- Indicator and "temperature" caption;

The back panel has:

➤ Nameplate with the indication of:

- Product name;
- Model;
- "In vitro diagnostics medical device" (the corresponding symbol is

marked **IVD**);

- Power supply conditions, input power;
- Product standard code;

- Production number;
- Production date;
- Manufacturer's name, address, telephone;
- Authorized official's name, address, telephone;
- Product license number and date,

Labels:

“Biohazard”

“Caution! Please, refer to the service manual”

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Fig.17 - ЮОНОН[®] LABSTAR 50 System main module nameplate marking

2.5.3 Transport (consumer) package marking:

- Manufacturer's name and address;
- Distributor's name and address;
- Product name, model;
- Serial number;
- Weigh net/ gross;
- Volume (length x width x height);
- Manufacture date;
- Product standard code;
- Indication “Fragile, handle with care” (the corresponding symbol is marked);
- Indication “Keep dry” (the corresponding symbol is marked);
- Indication “Up” (the corresponding symbol is marked);

- Indication “Do not roll” (the corresponding symbol is marked);
- Indication “Limit of the pile levels number”;
- Indication “Fragile precision tool. Caution”;
- transportation and storage conditions (the corresponding symbols are marked);

The designating of the marked symbols:

	–	Fragile, handle with care
	–	Keep dry
	–	Up
	–	Do not roll
	–	Limit of the pile levels number
	–	Temperature limit
	–	Humidity limit
	–	Atmosphere pressure limit

3 Medical device operation

3.1 Environmental requirements

3.1.1 Operational conditions

To provide the normal operation System should be installed and storage in the following environment:

- Environmental temperature: 10°C-30°C;
- Relative humidity: 45- 85%;
- Atmospheric pressure: 86.0 kPa ~106.0 kPa
- A/C power source: 220±22 V, frequency 50±1 Hz;

- Keep System away from intrusion of dust, bright directed light, liquids, corrosive gases, vibration and high electromagnetic disturbance

3.1.2 Transportation and storage requirements.

System is the accurate recording instrument. The special caution should be exercised to avoid the damage during transportation and storage. The device outer package should contain the corresponding marking labeling with the indications: "Handle with care", "Up", "Keep dry" etc. The carrier should take care to ensure that the device is not subjected to heavy shocks and to atmospheric precipitation or direct sunlight since the moment of shipment to delivering to the installation site. The following demands are placed on the transportation conditions:

- Environmental temperature: $-20^{\circ}\text{C} \sim +55^{\circ}\text{C}$;
- Relative humidity: 45~95%;
- Atmospheric pressure: 50.0 kPa ~ 106.0 kPa.

3.1.3 Power requirements

The analyzer is designed for the AC power supply with the voltage of $220\text{ V} \pm 22\text{ V}$, the frequency of $50\text{ Hz} \pm 1\text{ Hz}$; The device input power is 450 VA. There must be the regulated power supply between the device and the power system. The electric fuse characteristics are indicated on the back panel. When replacing the electric fuse, it is important to take the power cable characteristics into account.

 The ground wire of the three-conductor power cable should be directly connected to the safe ground network with the resistance of min 0.1 Ω to avoid the accidents.

3.2 Installation procedure

All the installation and implementation works of System must be carried out by the specialist assigned by the manufacturer's authorized official in the Russian Federation for this purpose. Additionally, the user should carry out the device acceptance.

3.2.1 Unpacking and acceptance

Open the package, carefully remove System with the accessories from the packing box. System packing list is inside of the box. After the unpacking, please check the contents of delivery according to the packing list. In case of the deficiency, address the distributor immediately.

3.2.2 Installation and implementation

The preoperational phase should be carried out in accordance with the

requirements listed in this section.

3.2.2.1 Power line

System should be connected to the separate power line with the ground fault protection. The distance between the device back and the wall should be of no less than 20 cm.


Fig. 18 - LABSTAR 50 installation space diagram

3.2.2.2 The switch, electrical outlet and electric fuse

The main module of System is equipped with the key-operated switch for 15A/30A 250 V, triphase electrical outlet for 10A, 250V and the electric fuse for 5A. The electrical outlet and the electric fuse holder are made as the single unit. All these elements are located on the right side of the back panel of the main module.


Fig.19- Back panel of LABSTAR 50 main module
 1-3 – Fans, 4 – Switch, 5 - power connector + electric fuse,
 6 - RS232 port for connecting the data cable

3.2.2.3 Main module installation

To install the main module, it is required only to connect the data cable, the power cord and the accessories. The device does not have any parts requiring the on-site dismantling. Carefully remove the main module from the package and install it on a stable lateral support or the table. Open the door, remove all the accessories that were fixed in its inside for the transportability. Make sure there are no foreign items inside of the device. Connect the main module to the computer COM port (RS232) using the data cable, connect the power cable.


Fig.20 - LABSTAR 50 connection diagram

⚠ The data cable in LABSTAR 50 model is made as the three-wired multiway cable. The cable terminal is connected to the main module communications port.

One of three wires is connected to the computer, the other two – with the PCI computer data collection card interface.

3.2.2.4 Computer installation

The computer installation should be carried out in the commonly-accepted way.


3.2.2.5 Commencement of work and operation

Before switching on the Blood Measure System main module, check the connection of the signal and data cables and the power cord. Make sure that the device total rated current does not exceed the power-supply source rated current. After that, switch on the power supply and start out the main module. The Blood Measure System main module will start operation in automatic mode.

Switch on the computer, double click the screen and enter the program. The Blood Measure System main module will start operation in automatic mode.

The registration process that is carried out by the manufacturer's authorized official is required for the initial system start up.

In order to do that, enter the organization name in "Department name" field and the System code indicated in the delivery set accompanying documentation, into "Device code" field of the "Logging in" pop-up window. Then obtain the license number from the manufacturer and enter it in the "License" field.


Fig.21 - Logging in

After pressing the “Determination” button the device will be registered and ready for use.

3.3 Setting of system parameters

3.3.1 Culture time

The system automatically sets the culture time for 120 hours. If there are special requirements, the time parameter on the screen should be modified when the bottles are put in and the time set on the “Set time” screen.


Fig. 22 - Bottles installation

3.3.2 Culture temperature

Culture temperature is constantly set as $35 \pm 1.5^\circ\text{C}$, which cannot be modified.

3.3.3 Setting the department name

Enter the instrument control procedure, click “System” and the interface of the dialog box appears. Click “Department” and then click “Add” to input the name of the specific department. Click “Delete” to delete the chosen department.


Fig.23 - LABSTAR 50 "System" screen

3.3.4 System reset

Enter the instrument control procedure, click "System" and the interface of the dialog box appears. Click "Cleared bottles". It will reset all the bottle positions".


Fig. 24 - LABSTAR 50 "System" screen

⚠ System reset can never be used except that serious faults of the instrument occur or power failure lasts for more than 1 hour. It cannot be used when the instrument is in normal operation, otherwise, all the data on test will be lost!!!

3.4 Operation procedure

3.4.1 Blood and other body fluids sample collection

Blood and other body fluids sample collection should be carried out according to the rules adopted in the laboratory, using the special bottles

Test method: Check the quality using visual inspection.

B1.1.5.1 Terminal

Requirement: The terminals should be recognized.

Test method: Check the quality using visual inspection.

B1.1.5.2 Terminals for connection to the mains supply

Not applicable

B1.1.6 Switches

Requirement: It should meet the requirement of 5.1.6 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.1.7 Equipment protected with double insulation and reinforced insulation

Requirement: It should meet the requirement of 5.1.7 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.1.8 The junction box

Not applicable

B1.2 Warning sign

Requirement: It should meet the requirement of 5.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.3 Durability of mark

Requirement: The marks should meet the requirements of standard provisions 5.1.2-5.2 and keep clearly distinguishable under the normal conditions.

Test method: According to the method of 5.3 in GB4793.1-2007.

B1.4 File

B1.4.1 Overview

Requirement: It should meet the requirement of 5.4.1 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.4.2 List of rated values

Requirement: It should meet the requirement of 5.4.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.4.3 Equipment installation

Requirement: It should meet the requirements of 5.4.3 a), b), c) in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.4.4 Operation of the equipment

Requirement: The following instructions should be included:

- a) Identification and use of the control buttons in the working conditions.
- b) The explanation for the symbols.
- c) The instruction for connection.

Test method: Check the quality using visual inspection.

B1.4.5 Maintenance of the equipment

Requirement: The following contents should be given in the specification:

- a) Preventive maintenance and checking by user.
- b) The parts which only be checked by manufacturer.
- c) The rating and feature of the fuse.

Test method: Check the quality using visual inspection.

B2 Electrical shock protection

B2.1 Overview

B2.1.1 Requirement

Requirement: The electrical shock protection should be kept in the normal or single failure condition of the device. The palpable parts of the device must not become dangerous electricity parts.

Test method: Determine whether the parts are palpable according to the rules of 6.2. Determine whether the test value reaches the limited value according to the rules of 6.2. Check the quality according to the method of 6.4-6.11.

B2.1.2 Exception

Not applicable

B2.2 Judging of the palpable parts

Requirement: It should meet the requirement of 6.2 in GB4793.1-2007.

Test method: According to the method of 6.2 in GB4793.1-2007.

B2.2.1 Checking

Requirement: It should meet the requirement of 6.2.1 in GB4793.1-2007.

Test method: According to the method of 6.2 in GB4793.1-2007.

B2.2.2 Openings above the charged parts

Not applicable

B2.2.3 Openings of the premodulation controls

Not applicable

B2.3 Allowable limit of the palpable parts

Requirement: It should meet the requirement of 6.3 in GB4793.1-2007.

Test method: Determine whether conform to the limited value according to the method of 6.5.1-6.5.4 in GB4793.1-2007.

B2.3.1 Value under the normal condition

Requirement: The voltage limitation are 33V (rms value) and 46.7V (peak value).

Test method: According to the method of 6.3.1 in GB4793.1-2007.

B2.3.2 Limits in the single failure condition

Requirement: The voltage limitation are 55V (rms value) and 78V (peak value).
The rms value of the sine-wave current is 3.5mA.

Test method: According to the method of 6.3.2 in GB4793.1-2007.

B2.4 Protection in normal condition

Requirement: The following one or more measures should be applied to the palpable parts in order to prevent dangerous electricity:

—Basically insulation

—Shell or baffles

Test method: Check the quality using visual inspection.

B2.5 Protection in single failure condition

Requirement: The additional protection should be provided to ensure that the palpable parts cannot become dangerous electricity in single failure condition. The additional protection is composed of one or more measures specified by 6.5.1-6.5.3.

Test method: Determine whether conform to the rules according to the methods of 6.5.1-6.5.4 in GB4793.1-2007.

B2.5.1 Protective connection

Requirement: If the primary protection equipment specified by 6.4 could lead to conducting of the electric parts in single failure condition, then the palpable parts should be connected with protective conductor terminal.

Test method: Determine whether conform to the rules according to the method of 6.5.1 in GB4793.1-2007.

B2.5.1.1 Integrity of the protective connection

Requirement: It should meet the requirement of 6.5.1.1 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.5.1.2 Protective conductor terminal

Requirement: It should meet the requirement of 6.5.1.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection and according to the method of 6.5.1.3 in GB4793.1-2007.

B2.5.1.3 Connection impedance of the plug connection device

Requirement: The impedance of the protective connection which is used for protective conductor terminal between each accessible part shall not exceed 0.1Ω .

Test method: According to the method of 6.5.1.3 in GB4793.1-2007.

B2.5.1.4 Connection impedance of the permanently-connected device

Not applicable

B2.5.1.5 Screen as the transformer protective bonding

Not applicable

B2.5.2 Double insulation and reinforced insulation

Requirement: It should meet the requirement of 6.5.2 in GB4793.1-2007.

Test method: According to the method of 6.5.2 in GB4793.1-2007.

B2.5.3 Protective impedance

Not applicable

B2.5.4 Automatic disconnecting of the power supply

Not applicable

B2.6 Connection with an external circuit

B2.6.1 Overview

Requirement: It should meet the requirement of 6.6.1 in GB4793.1-2007.

Test method: According to the method of 6.6.1 in GB4793.1-2007.

B2.6.2 External circuit terminal

Not applicable

B2.6.3 Circuit with dangerous charged terminal

Not applicable

B2.6.4 Palpable terminal for stranded conductor

Not applicable

B2.7 Requirements to insulation: Clearance and creep age distance

Creep age distance: insulation type: basic insulation and reinforced insulation.

Microenvironment pollution level: grade 2

Requirement: 1.5mm

Creep age distance: insulation type: basic insulation. II Material group:

$400 \leq C_{IT} \leq 600$.

Requirement: 3.0mm

Creep age distance: insulation type: basic insulation. II Material group:

$400 \leq C_{IT} \leq 600$.

Requirement: 6.0mm

Test method: Verified by checking

B2.8 Dielectric strength test

Requirement:

1) The basic insulation part: It should be tolerated the test in the condition of 50Hz and 1390V. The breakdown or repeating flash-over must not appear.

2) The double insulation or reinforced insulation part: It should be tolerated the test in the condition of 50Hz and 2224V. The breakdown or repeating flash-over must not appear.

Test method: The voltage test is measured according to the specified voltage value of 6.8.4 in GB4793.1-2007. The breakdown or repeating flash-over must not appear. The corona discharge effect or similar phenomenon is negligible.

B2.9 Structure requirements of protection against electric shock

B2.9.1 Overview

Requirement: It should meet the requirement of 6.9.1 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.9.2 Insulation materials

Requirement: It should meet the requirement of 6.9.2 in GB4793.1-2007.

Test method: According to the method of 6.9.2 in GB4793.1-2007.

B2.9.3 Color coding

Not applicable

B2.10 Connection to electric mains and between different parts

B2.10.1 Power cords

Requirement: It should meet the requirement of 6.10.1 in GB4793.1-2007.

Test method: According to the method of 6.10.1 in GB4793.1-2007.

B2.10.2 Installation of the unable-dismantled power cords

Not applicable

B2.10.3 Plug and connector

Requirement: It should meet the requirement of 6.10.3 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.11 Disconnection of the power supply

B2.11.1 Overview

Requirement: The device should be equipped with breaking device which can turn off the power supplied from each electric source either inside or outside of the device. The breaking device should disconnect all the current carrying conductors.

Test method: According to the method of 6.11.1.1-6.11.3.2 in GB4793.1-2007.

B2.11.1.1 Exception

Not applicable

B2.11.2 Requirements according to the type of the device

B2.11.2.1 Permanently connection and multi-phase equipment

Not applicable

B2.11.2.2 Single phase equipment connected with the cord

Requirement: It should meet the requirement of 6.11.2.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.11.2.3 Dangerous caused by functions

Not applicable

B2.11.3 Breaking device

Requirement: It should meet the requirement of 6.11.3 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.11.3.1 Switch and circuit breaker

Requirement: It should meet the requirement of 6.11.3.1 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B2.11.3.2 Appliance coupler and plug

Requirement: It should be meet the requirement of 6.11.3.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B3 Prevention of mechanical hazard

B3.1 Overview

Requirement: It should meet the requirement of 7.1 in GB4793.1-2007.

Test method: According to the method of 7.2-7.6 in GB4793.1-2007.

B3.2 Moving machine part

Requirement: It should meet the requirement of 7.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B3.3 Stability

Requirement: It should meet the requirement of 7.3 in GB4793.1-2007.

Test method: According to the method of 7.3 in GB4793.1-2007.

B3.4 Lift and handling device

Requirement: It should meet the requirement of 7.4 in GB4793.1-2007.

Test method: According to the method of 7.4 in GB4793.1-2007.

B3.5 Installation on the wall

Not applicable

B3.6 Expelled parts

Not applicable

B4 Resistance to mechanical impact

Requirement: After the test is complete, the equipment should pass the voltage test according to 6.8.4, and by visual inspection to verify:

- The dangerous charged components should not become accessible.
- There should not have any dangerous cracks on the shell.
- The clearance should not less than the allowable values, and the insulation of the internal wire should not be damaged.
- The baffle should not be damaged or loose.
- Motion components cannot be exposed except those allowed in 7.2.

-- There should be no damage caused by the spread of flame.

Test method: According to the method of 8.1 and 8.2 in GB4793.1-2007.

B4.1 The rigid test of the shell

B4.1.1 Static test

Requirement: It should meet the requirement of 8.1.1 in GB4793.1-2007.

Test method: According to the method of 8.1.1 in GB4793.1-2007.

B4.1.2 Dynamic test

Requirement: It should meet the requirement of 8.1.2 in GB4793.1-2007.

Test method: According to the method of 8.1.2 in GB4793.1-2007.

B4.2 Drop test

Not applicable

B5 Prevent the spread of flame

Under normal conditions or a single fault conditions, the flame shall not be spread to the outside of the equipment. Test whether it is qualified or not by at least one of the following methods.

a) Do the test under a single fault condition which could lead to the spread of flame to the equipment outside.

b) Check whether the flame can be controlled inside the device once a fire happens according to the provisions in 9.2.

Test method: Verified by checking

B5.1 Eliminate or reduce ignition sources within the equipment

Not applicable

B5.2 Control the flame inside the equipment once a fire happens.

The equipment and its shell should meet the structure requirements in 9.2.1.

Test method: According to the method of 9.1 in GB4793.1-2007.

B5.3 Energy-limited circuits

Not applicable

B5.4 Requirements of the device containing or using flammable liquids

Not applicable

B5.5 Over current protection

Predetermine that the device powered by electric network source should be protected by the using of fuse, circuit breaker, thermal circuit breaker and impedance limiting circuit. OCP device should also provide protection against electric shock under fault conditions.

OCP device should not be installed in the protective wire. The fuse or single-pole circuit breaker should not be installed in the middle wire of the multi-phase device.

Test method: According to the method of 9.5 in GB4793.1-2007.

B5.5.1 Permanent connection device

Not applicable

B5.5.2 Other device

Requirement: The OCP (over current protection) device should be installed inside the equipment if it is adopted.

Test method: Check the quality using visual inspection.

B6 Temperature limits and heat-resisting

B6.1 Surface temperature limit of burning prevention

Not applicable

B6.2 Winding temperature

Not applicable

B6.3 Measurement of other temperature

Requirement: At 40 degrees or maximum rated of the ambient temperature (if more higher), the clearance and creep age distance should meet the requirements of 6.7.

Test method: According to the method of 10.3 in GB4793.1-2007.

B6.4 Implementation of temperature test

Test method: According to the method of 10.4 in GB4793.1-2007.

B6.4.1 Temperature test of heating element

Not applicable

B6.4.2 Device designed for installation in the cabinet or on the wall

Not applicable

B6.5 Heat resistance

B6.5.1 Integrity of clearance and creep age distance

Requirement: It should meet the requirements of 10.5.1 in GB4793.1-2007.

Test method: According to the method of 10.5.1 in GB4793.1-2007.

B6.5.2 Non-metallic case

Requirement: It should meet the requirement of 10.5.2 in GB4793.1-2007.

Test method: According to the method of 10.5.2 in GB4793.1-2007.

B6.5.3 Insulation material

Requirement: It should meet the requirement of 10.5.3 in GB4793.1-2007.

Test method: According to the method of 10.5.3 in GB4793.1-2007.

B7 Prevention of liquid dangerous

B7.1 Overview

Requirement: It should meet the requirements of 11.1 in GB4793.1-2007.

Test method: Verified by checking

B7.2 Cleaning

Requirement: It should meet the requirement of 11.1 in GB4793.1-2007.

Test method: Verified by checking

B7.3 Splashing

Requirement: It should meet the requirements of 11.2 in GB4793.1-2007.

Test method: Verified by checking

B7.4 Overflowing

Not applicable

B7.5 Battery electrolyte

Not applicable

B7.6 Special protective device

Not applicable

B7.7 Liquid pressure and leakage

Not applicable

B8 Protection from radiation (including the laser source), acoustic pressure and ultrasonic pressure

Not applicable

B9 Protection from the release of the gases and substances, explosion and implosion

Not applicable

B10 Components and assemblies

B10.1 Overview

If the security is involved, the component should be used according to the rated value unless there has been made certain exceptions.

Test method: Check the quality using visual inspection.

B10.2 Electric motor

B10.2.1 Electric motor temperature

Not applicable

B10.2.2 Series-excited motor

Not applicable

B10.3 Over-temperature protection device

Not applicable

B10.4 Fuse holder

Requirement: It should meet the requirement of 14.4 in GB4793.1-2007.

Test method: According to the method of 14.4 in GB4793.1-2007.

B10.5 Selective device of power supply voltage

Not applicable

B10.6 High sophistication component

Not applicable

B10.7 Power transformer testing outside the device

Not applicable

B10.8 Printed wiring board

Requirement: FV-1 (flammability rating GB/T 11020) or much better material should be used.

Test method: According to the method of 14.8 in GB4793.1-2007.

B10.9 Electrical circuit and component used for transient overvoltage limiting device

Not applicable

B11 Protection with interlock device

Not applicable

B12 Hazards resulting from application

Not applicable

10 Appendix C

Safety requirement, test method and inspection rule (according to the YY 0648-2008 standard)

C0 Application condition of the instrument

Application condition of the instrument

- a) For indoor use.
- b) Temperature: 10-30°C
- c) Humidity: No more than 85%
- d) Power condition: The power fluctuation is no more than 10% of the standard voltage
- e) Overvoltage category: The transient overvoltage category II
- f) Pollution level: Grade II

C1 Mark and file

C1.1 Mark

C1.1.1 Overview

Requirement: The words and symbols marked on the outside should meet the requirements of 5.1.1 (Table 1) in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.1.2 Mark

Requirement: It should meet the requirements of 5.1.2 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.1.3 Transportation and storage

Requirement: It should meet the requirements of 5.1.101 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.2 Warning sign

Requirement: It should meet the requirements of 5.2 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.3 Durability of the marks

Requirement: Marks required in 5.1.2-5.2 of YY 0648-2008 should be permanently attached and keep clear under the condition of normal use.

Test method: Check according to the method in 5.3 of YY 0648-2008.

C1.4 File

C1.4.1 Overview

Requirement: It should meet the requirements of 5.4.1 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.4.2 Transportation, installation and assembly instruction of the instrument

Requirement: It should meet the requirements of 5.4.3 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.4.3 Operation of the instrument

Requirement: It should meet the requirements of 5.4.4 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.4.4 Stop using of the device because of maintenance or handling

Requirement: It should meet the requirements of 5.4.101 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C1.4.5 Transportation and storage

Requirement: It should meet the requirements of 5.4.102 in YY 0648-2008.

Test method: Check the quality using visual inspection.

C2 Transportation and storage

Requirement: It should meet the requirement of 8.101 in YY 0648-2008.

Test method: According to the method of 8.101 in YY 0648-2008.

C3 Protection for the release of the gas, explosion and implosion

Not applicable.

11 Appendix D

Automated bacteriological system **ЮHOHA®** LABSTAR, design versions:

ЮHOHA® LABSTAR 50, **ЮHOHA®** LABSTAR 100 with accessories meets the requirements for emissions and noise immunity GB/T 18268.1-2010 and GB/T 18268.26-2010.

The user must ensure that the analyzer is used in a suitable electromagnetic environment.

Use of other electrical equipment on or near the analyzer may cause interference. Before using the analyzer for its intended purpose, you must verify that the configuration described in this Guide is correct.

Do not use the analyzer near strong radiation sources (such as unshielded RF sources) as this may interfere with normal operation of the device.

Automated bacteriological system **ЮHOHA®** LABSTAR, design versions:

ЮHOHA® LABSTAR 50, **ЮHOHA®** LABSTAR 100 with accessories designed for use in basic electromagnetic environment as below.

Electromagnetic emission

Electromagnetic emission test	Conformity	Electromagnetic environment - instructions
Radio emission	Group 1	Automated bacteriological system ЮHOHA® LABSTAR, design versions: ЮHOHA® LABSTAR 50, ЮHOHA® LABSTAR 100, with accessories uses radio emission only to perform its internal function. Thus, the radio emission of the product is extremely small and probably cannot cause interference electronic equipment located nearby.
Radio emission	Class A	Automated bacteriological system ЮHOHA® LABSTAR, design versions: ЮHOHA® LABSTAR 50,

		ЮHOHa® LABSTAR 100, with accessories suitable for use in all locations except residential buildings and buildings directly connected to the distribution network supplying residential buildings. When used in a residential area, it may cause radio interference. In such cases, the user should take appropriate measures.
Harmonic currents	Class A	
Voltage fluctuations and flicker	Corresponds	

Electromagnetic interference

Port	Name of electromagnetic interference	Standard EMC	Test parameter value	Criterion of quality of functioning
The port of hull	Electrostatic discharge (ESD)	GB/T 17626.2 (IEC 61000-4-2)	±4 kV (contact discharge) ±8 kV (air discharge)	A A
	Electromagnetic field	GB/T 17626.3 (IEC 61000-4-3)	3 V / m (80 MHz to 1 GHz) 3 V / m (1.4 GHz to 2 GHz) 3 V / m (2 GHz to 2.7 GHz)	A
	Industrial frequency electromagnetic field	GB/T 17626.4 IEC 61000-4-8	3 A / m (50 Hz)	A

Port AC power (including protective earth ports)	Voltage sags	GB/T 17626.11 (IEC 61000- 4-11)	0%, first half of the period	A
			0%, 1 period	A
			70%, 25 periods	A
	Short-term power supply voltage interruptions	GB/T 17626.11 (IEC 61000- 4-11)	0%, 250 periods	C
	Nanosecond impulse noise	GB/T 17626.4 (IEC 61000- 4-4)	1 kV (5/50 NS, 5 kHz)	A
	Microsecond pulses of high energy	GB/T 17626.5 (IEC 61000- 4-5)	1 kV Line-ground	A
	Conductive interference induced by radio frequency fields	GB/T 17626.6 (IEC 61000- 4-6)	3 V (150 kHz to 80 MHz)	A

Criterion of quality of functioning A - the Equipment must continue to function as intended during and after the test. Deterioration of performance or loss of functions below the level of quality of operation provided in the technical documentation of the manufacturer for the equipment when using the equipment for its intended purpose is not allowed.

Criterion of quality of functioning C - temporary loss of functions is allowed provided that the function can be restored independently or during operation of controls.

12 Appendix E

List of applicable national and international standards to ensure the safety, effectiveness and quality of the medical device

The medical device "Automated bacteriological system **ЮHOHA®** LABSTAR, design versions: **ЮHOHA®** LABSTAR 50, **ЮHOHA®** LABSTAR 100 with the accessories" corresponds to the national standards of the People's Republic of China for the products:

- GB 4793.1-2007 "Safety requirements for electrical equipment for measurements, control and use in laboratories - Part 1: General requirements";
- GB/T 14710-2009 "Environmental requirements and the method for controlling medical electrical equipment";
- YY 0466-2003 "Symbols of medical devices used for labels and markings of medical devices or for the provision of information";
- YY/T 0648-2008 "Safety requirements to electrical equipment for measurement, control and use in laboratories - Part 2-101: Specific requirements to medical devices used for in-vitro diagnostics (IVD)";
- YY/T0316-2008 "Medical devices. Application of risk management for medical devices";
- YY/T 0656-2008 "Automated Blood Culture System";
- GB/T 191-2008 "Packaging - Marking of transport cargoes";
- GB/T 25000.51-2010 "Requirements to software development and quality, testing methods".
- GB/T 18268.26-2010 «Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 26: Particular requirements. In vitro diagnostic (IVD) medical equipment»;
- GB/T 18268.1-2010 «Electrical Equipment for Measurement, Control and Laboratory Use - EMC Requirements - Part 1: General Requirements»;
- GB 4824-2013 (CISPR 11: 2003) «Industrial, scientific and medical (ISM) radio-frequency equipment. Disturbance characteristics. Limits and methods of measurement»;
- GB/T 17626.3-2016 «Electromagnetic compatibility- Testing and measurement techniques- Radiated , radio-frequency , electromagnetic field immunity test»;

- GB/T 17626.4-2018/IEC 61000-4-4:2012 «Electromagnetic compatibility. Testing and measurement techniques. Electrical fast transient/burst immunity test»;
- 17626.5-2008/IEC 61000-4-5:2005 «Electromagnetic compatibility. Testing and measurement techniques. Surge immunity test»;
- GB/T 17626.6-2017 «Electromagnetic compatibility. Testing and measurement techniques. Immunity to conducted disturbances, induced by radio-frequency fields»;
- GB/T 17626.8-2006 «Electromagnetic compatibility (EMC). Part 8: Testing and measurement techniques. Power frequency magnetic field immunity test»;
- GB/T 17626.11-2008 «Electromagnetic compatibility. Testing and measurement techniques. Voltage dips, short interruptions and voltage variations immunity tests».

The medical device “Automated bacteriological system **ЮХОХа®** LABSTAR, design versions: **ЮХОХа®** LABSTAR 50, **ЮХОХа®** LABSTAR 100 with the accessories” corresponds to the national standards of the Russian Federation for the products:

- GOST R 50444-92 «Medical instruments, apparatus and equipment. General specifications»;
- GOST 12.2.091-2012 (IEC 61010-1:2001) «Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1. General requirements»;
- GOST IEC 61010-2-081-2013 «Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 2-081. Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes»;
- GOST IEC 61010-2-101-2013 «Safety requirements for electrical equipment for measurement, control and laboratory use. Part 2-101: Particular requirements for in vitro diagnostic (IVD) medical equipment»;
- GOST IEC 60825-1-2013 «Safety of laser products. Part 1. Equipment classification, requirements and user's guide»;

- GOST R ISO 18113-1-2015 «In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements»;
- GOST R ISO 18113-3-2015 «In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 3. In vitro diagnostic instruments for professional use»;
- GOST R EN 13612-2010 «Performance evaluation of in vitro diagnostic medical devices»;
- GOST R IEC 62304-2013 «Medical devices. Software. Life cycle processes»;
- GOST R ISO/IEC 12119-2000 «Information technology. Software packages. Quality requirements and testing»;
- GOST R ISO/IEC 9126-93 «Information technology. Software product evaluation. Quality characteristics and guidelines for their use»;
- GOST 28195-89 «Quality control of software systems. General principles»;
- GOST R IEC 62366-2013 «Medical devices. Application of usability engineering to medical devices»;
- GOST R 56894-2016/GHTF/SG1NO63:2011 «Summary technical documentation for demonstrating conformity to the essential principles of safety and performance of in vitro diagnostic medical devices»;
- GOST R IEC 61326-1-2014 «Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1. General requirements»;
- SanPiN 2.1.7.2790-10 «Sanitary and epidemiological requirements for the treatment of medical waste».

manufactured by SCENKER Biological Technology Co., Ltd. (see the more detailed information about the analysis of blood and other body fluids in the Cultural medium Bottles documentation).

- The blood and other body fluids sample for the culture should be taken with keeping the rules of the aseptic to avoid the microorganisms' penetration from the outside environment.

- Before using the bottles, medium transparency should be visually assessed, any opacity means their unfitness.

- The bottles should be warmed for the room temperature (18-30°C).

- The culture should be carried out during the temperature rise, when the fever begins. To take blood using special sterile systems for blood sampling or sterile syringes in the bottle created a vacuum, which facilitates the sowing of blood samples in the bottle. Break off the plastic covers from the bottles and sterilize the inside cover with the drape saturated in 70% ethanol (or other substitutive antiseptic) within one minute. Carefully sterilize the skin over the needled vein with 70% ethanol and then with 1-2% tincture of iodine within 30 seconds. After drying of the treated area, perform venipuncture. For adult patients it is recommended to use a pair of aerobic and anaerobic bottle. When using the system for blood sampling, the inoculation is performed first in an aerobic vial and then in an anaerobic one. (Do not remove the needle from the patient's vein during blood collection. After the end of the first extract needle from vial, then from veins patient). When using a syringe, first inject blood into an anaerobic vial and then into an aerobic one (Avoid air entering the anaerobic bottle). The volume of the blood is from 5 ml to 10 ml for the adults for aerobic bottle and from 2 ml to 10 ml for the anaerobic and children bottle (the marking on the name plate of the bottle allows you to estimate the volume of the taken blood). Optimal pathogen isolation is achieved by adding the maximum sample volume Mix the bottle contents carefully with the circular motions. After finishing the procedure, the iodine should be removed from the skin with the drape using the 70% ethyl and the puncture area should be covered the adhesive plaster. The bottles should be marked with indicating the name, initials, date of birth and sampling date, and then transported to the laboratory in the at the room temperature (18-30°C) in a dark place (not in the freezer) for no more than 24 hours.

- Cerebrospinal fluid is obtained with lumbar puncture, less often with puncture of the lateral ventricles of the brain, preferably before taking antibiotics. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of

iodine within 30 seconds. The first portions of liquor have the highest potential for contamination by microorganisms from the skin and should be discarded. Immediately after the puncture, not less than 5 ml and not more than 10 ml of the resulting liquor should be aseptically introduced into an aerobic vial (in certain clinical situations, an anaerobic vial is required). Large volumes (5-10 ml) increase the sensitivity of the culture and are necessary for optimal isolation of whimsical bacteria and fungi. If the volume of liquor is small from 2 to 5 ml, as well as liquor from children - you should use aerobic bottles for children. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of pathogens. During storage and transportation, excessive overheating and especially cooling should be avoided (a thermos or a thermal container with a temperature of $35 \pm 1.5^{\circ}\text{C}$ can be used). Cooling the liquor below 30°C leads to the death of meningococci. Storage and transportation time should not be more than 8 hours.

- In the study of body fluids from normally sterile cavities (peritoneal, synovial, pleural, articular, pericardial fluid) for adults, it is recommended to take for each aerobic vial from 5 ml to 10 ml, for each anaerobic vial from 2 ml to 10 ml of the test liquid, for children for each vial it is necessary from 2 ml to 10 ml of the test liquid. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. Immediately after the puncture, the biological fluid should be aseptically introduced into the aerobic vial (in certain clinical situations, an anaerobic vial is required). The minimum volume of liquid for bacteria isolation is 2-5 ml, for fungi isolation-10 ml. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of microorganisms. Store at room temperature ($18-30^{\circ}\text{C}$) for up to 24 hours.

- The rate of growth of pathogens is directly proportional to the volume of biomaterial. Optimal pathogen detection is achieved with maximum sample volume.

- After sampling and inoculation of blood or biological fluids, it is necessary to tear off half of the label with a bar code from the bottle (the other half of the label can not be torn off) and stick it on the form of the referral for laboratory research

When separating the label with the barcode, it is necessary to check that the information on the issue is the same as on the bottle!

⚠ The samples should be forwarded to the laboratory for testing as soon as possible. This will increase the possibility of the pathogenic microorganism detection.

3.4.2 Bottles loading

3.4.2.1 Using the barcode scanner

Click “Open” in the control procedure and open the door of the main module. After the door is open, “Open” will be inactive.



Fig. 25 - Opening the door before the bottle installation

Click “Load”. The dialog box will automatically pop up on the screen. Use the barcode scanner to read the barcode from the label. System will automatically allocate the Position No. for placing the bottle. Put the bottle to the necessary position and click “Load” (at the same time the Position No. at main program interface manifests yellow) (fig.24). Load all the bottles for the check and click “Confirm”. Close the main module door and it will automatically resume operation.


Fig. 26 - Loading the bottle

3.4.2.2 With the manual entering the barcode

Open the door of the main module, click “Load”. The dialog box will automatically pop up on the screen. Manually enter the barcode from the label. System will automatically allocate the Position No. for placing the bottle. Put the bottle to the necessary position and click “Load” (at the same time the Position No. at main program interface manifests yellow). Load all the bottles for the check and click “Confirm”. Close the main module door and it will automatically resume operation.

⚠ During the operation of the main module it is not recommended to often open the door or hold it open. It can cause the temperature fluctuations inside of the module and affect the test results

3.4.3 Input, inquiry and statistics of patients operation

Information Input: Enter the control procedure of the instrument, click “Data”, and the interface of dialog box will appear. Click “Specimens Information” and use the bar code scanner to read the bar code. The instrument will generate the ID No. (It is the same as that of Bar Code No.) Following that, fill the patients data (Surname, first name & patronymic, gender, age, department name, patient number, patient history No., date of check) in the dialog box. Click “Save” after patient data is filled.


Fig.28 - Information inquiry


Fig. 29 - Information statistics by date

Besides, the detection results can be displayed in the graphic view.


Fig.30 - Information inquiry, growth curve

The graphic image can be saved on the computer: Click "Printing plan" and the

dialog box for saving the document will appear


Fig. 31 - Saving the document

The image is saved in mdi format.

3.4.4 Culture

After the culture bottle is put into the main module, System will automatically start the culture. The temperature should be maintained at the level of $35 \pm 1.5^\circ\text{C}$. The LABSTAR 50 turntable will rotate with the constant speed of 26 r/min.

3.4.5 Detection

The instrument operation cycle for the dynamic detection lasts for 10 minutes. The culture bottle stimulates the generated fluorescence as affected by the bacteria metabolism products, and sends the corresponding signals to the analysis program of the system computer. The provision is made for the previous transformation of the analog signal into the digital one. The Position No. of the culture bottle in the main program interface manifests yellow.

3.4.6 Alerting system

The system will analyze the detection signal. It will alarm immediately when it discovers the positive bottle. The Position No. of the culture medium in the main program interface manifests red after the alarm. The figure shows the alarming point.


Fig. 32 - Growth curve with alarming point for the positive bottle

Unload the positive bottle from the device and use the sample for the subculture to conduct the additional analysis experiments. The culture bottle with no growth of culture medium of microorganism after 5 days (120 hours) of successive culture will be reported as the negative bottle and the corresponding alarm signal will occur. The Position No. of the culture bottle in the main program interface will manifest grass-green and the culture bottle can be taken out.

⚠ After all the negative bottles are taken out, the sub-culture should be applied immediately and the additional analysis should be carried out so as to prevent the instrument from alarming falsely!

3.4.7 Unloading the bottles

3.4.7.1 Using the barcode scanner:

Open the device door and unload the negative bottles (green color on the interface screen) or the positive bottles (red color on the interface screen). Click “Unload”, the instrument will automatically pop up the dialog box, then use the bar code scanner to read the sample bar code. In the meanwhile, the Position No. of the culture bottle in the main program interface manifests white, which indicates that the culture bottle has been taken out.


Fig. 33 - Unloading the bottle

In case of the wrong or the early bottle unloading the system will pop up the alert message:


Fig.34 - Alert message

3.4.7.2 Manual entering the bar code:

Open the device door and unload the bottles from the corresponding position numbers (green or red colors on the screen). Click “Unload”, the instrument will automatically pop up the dialog box, then enter the bar code manually and click “Unload”. In the meanwhile, the Position No. of the culture bottle in the main program interface manifests white, which indicates that the culture bottle has been taken out.

⚠ After the culture bottle is taken out, please immediately close the instrument door, so as to prevent low temperature inside the instrument!

After finishing the operation at clicking “Close” the program is closed but System works on.


Fig. 35 - Closing the program

The short System workflow is shown on Figure 36.


Fig. 36 - LABSTAR System operation cycle flow diagram

3.5 Culture media and application scope

3.5.1 Culture media

Table 9

Name	The volume of the collected material	The volume of the medium	Application	Barcode/ ID No.	Photo
Culture medium Bottles ЮHOHA® for aerobe cultivation;	5-10 ml	25 ml	The samples of the patients who did not take the antibiotics	XBZ	
Antibiotic-neutralized Culture medium Bottles ЮHOHA® for aerobe cultivation;	5-10 ml	25 ml	The samples of the patients who took the antibiotics	XCR	
Antibiotic-neutralized Culture medium Bottles ЮHOHA® for anaerobe cultivation	2-10 ml	25 ml	The samples of the patients who took the antibiotics	XYX	
Antibiotic-neutralized Culture medium Bottles ЮHOHA® for aerobe cultivation for children.	2-10 ml	25 ml	The samples of the children or other patients that are difficult to collect	XET	

3.5.2 Application scope

The medium bottles are designed for the culture and detection of pathogens in the blood samples and other aseptic body fluids.

3.6 How to collect the blood and other sterile biological fluid samples

3.6.1 Procedure for blood sampling

3.6.1.1 Blood sample collection indication

As a rule, the indication for the corresponding clinical examination is the emergence of one or few symptoms of the following list:

- Fever ($\geq 38^{\circ}\text{C}$) or low body temperature ($\leq 36^{\circ}\text{C}$);
- Shivering;
- Leukocyte increasing (more than $10.0 \times 10^9/\text{l}$, especially for the “left deviation”)
- Mucocutaneous hemorrhage;
- Coma;
- Multiple organ failure;
- Blood pressure decreasing;
- Increased C-reactive protein and breath accelerating;
- Blood disease patients have the reduction of granulocyte and platelet, etc.

3.6.1.2 Examination frequency

- The collection usually takes place at the early or peak fever periods
- The patient who is suspected to have the bacterial endocarditis should pass three examinations to determine the pathogens and the blood collection is carried out no less than one time per hour within 24h.

➤ If the patient was treated with antibiotics in one to two weeks before the culture, it is important to carry out two or three microbiological studies within one day afterwards.

➤ Every time the sample should be divided into two or three parts. The collection is carried out within 24h.

3.6.1.3 Collected material volume

The collected material volume corresponds to the Table 9

3.6.1.4 Skin sterilization process

- Sterilize the skin over the needled vein with 70% alcohol solution within 30 seconds;
- Sterilize with 1-2% tincture of iodine or tamed iodine within 60 seconds;
- Remove the rests of iodine from the skin with 70% alcohol solution.
- Carry out the blood collecting after the alcohol evaporation.

3.6.1.5 Bottle sterilization process

- Sterilize the culture bottle rubber stopper with 70% alcohol solution within 60 seconds;

➤ Remove the alcohol rests from the rubber stopper with the surgical gauze or cotton stopper.

➤ Insert the blood samples to the bottle.

3.6.1.6 Blood samples storage and transportation

➤ The samples should be sent to the laboratory for the examination after the blood collection;

➤ If it is impossible to send the sample for the examination at once, it can be stored at the room temperature (18-30°C) or in the constant-temperature bath (35±1,5°C). Cold storage is not allowed;

➤ The delay for the examination cannot exceed 24 h after collection.

3.6.1.7 Laboratory receiving standards

➤ Check the samples for the visible turbidity, pellicle, bacterial lawn, coagulation and other bacterial development signs.

➤ Check whether the blood volume is suitable.

➤ Check the bottle for the leaking, cracks and visible contaminations.

➤ Check whether the label on the bottle is in compliance with the application list and whether the patients' information is complete.

3.6.1.8 Samples screening procedure

➤ Refuse to take the bottles for the examination if the information on the bottle is not in compliance with the assay application and the culture bottles are broken or have the visible contaminations.

➤ The bottles are rejected if they are kept for more than 24 h after the collection, the bottle type does not correspond to the examination type, the bottle is expired, the volume is not enough for the examination etc.

3.6.2. Cerebrospinal fluid is obtained with lumbar puncture, less often with puncture of the lateral ventricles of the brain, preferably before taking antibiotics. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. The first portions of liquor have the highest potential for contamination by microorganisms from the skin and should be discarded. Immediately after the puncture, not less than 5 ml and not more than 10 ml of the resulting liquor should be aseptically introduced into an aerobic vial (in certain clinical situations, an anaerobic vial is required). Large volumes (5-10 ml) increase the sensitivity of the culture and are necessary for optimal isolation of whimsical bacteria and fungi. If the volume of liquor is small from 2 to 5 ml, as well as liquor from children - you should use aerobic bottles for children. Samples should be sent to the laboratory for examination as soon as possible. This will increase the

probability of detection of pathogens. During storage and transportation, excessive overheating and especially cooling should be avoided (a thermos or a thermal container with a temperature of $35 \pm 1.5^{\circ}\text{C}$ can be used). Cooling the liquor below 30°C leads to the death of meningococci. Storage and transportation time should not be more than 8 hours.

3.6.3. In the study of body fluids from sterile normal cavities (peritoneal, pleural, articular, pericardial fluid) for adults is recommended to take for each aerobic vials from 5 ml to 10 ml, for each anaerobic vial of 2 ml to 10 ml of the test liquid, for children of each vial must be from 2 ml to 10 ml. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. Immediately after the puncture, the biological fluid should be aseptically introduced into the aerobic vial (in certain clinical situations, an anaerobic vial is required). The minimum volume of liquid for bacteria isolation is 2-5 ml, for fungi isolation-10 ml. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of microorganisms. Store at room temperature ($18-30^{\circ}\text{C}$) for up to 24 hours.

3.7 Results explanation

The human blood and body fluids usually has no bacteria. The bacterial invasion can cause the bacteremia or septicemia. The patients with endocarditis, endarteritis, early typhoid and brucellosis are considered to be suffering from the bacteremia. Some other infections may be accompanied with intermittent bacteremia. Various results of the blood culture should be analyzed according to the specific situation. The culture results do not provide the exact bacteremia indication.

3.7.1 Negative results

The negative results can be explained in the following ways:

(1) There is no bacteremia

(2) The patient has the clinical signs but the bacteremia was not detected because of the following reasons:

- Local infection
- Fastidious bacteria, chlamydia and L-type bacteria samples
- Incorrect collecting time
- Using the antibiotic therapy before the examination
- Periodic bacteremia

Treatment measures (2):

- Blood samples re-collection and examination.

3.7.2 Positive results:

(1) Once positive:

- Pathogenic bacteria: Brucella, Salmonella ect. Bacteremia
- Saprophytes: Propionibacterium, Staphylococcus epidermidis

Contamination

(2) Several times positive:

- The same bacterium: meningitis, endocarditis and bacteremia
- Different bacterium (due to immunosuppression): dermatitis and bacteremia

4 Maintenance and repairs

4.1 The user-performed technical maintenance:

Type of the maintenance	LABSTAR 50
Cleaning the air filters with the fresh water	Once per 2 weeks
Dust collection in System functional area	Once per 2 months
Cleaning the liquids spilled on System	In case of the leakage or the bottle damage, remove the spilled liquid according to the decontamination procedure recommended in the laboratory. After the disinfection, clean with the wet wipe (moistened with water only) and rub dry.
Cleaning the liquids spilled inside of System	Remove the samples spilled inside of System and carry out the disinfection according to the authorized procedure.
Checking System cleanness	Once per 3 months, clean with the dry cotton plug if necessary
Checking the output voltage of the regulated power supply	Once per 6 months
Fixed temperature and humidity maintenance in the laboratory	Within 24 hours
System supply disconnection in case of the power outage	Turn on System after the power restart
Switch off System power supply and restart it in case of the sputter of the signal system	In 3 minutes after shutdown

and failure to correct the error

4.2 The common failures and treatment measures

Error category	Cause	Solutions
Main module power fail	No power cable or the power cable unplugged.	Check the power cable connection at both ends.
	The main module power switch is in the "off" position.	Set the power switch into "On" position (it is located in the lower right corner on the main module back panel).
Lost connection between the PC and the main module	No data cable between the PC and the main module or it is wrongly connected.	Check the cable connection at both ends. Restart the main module.
Temperature Alarm	The temperature of the incubator is too high	1. Reduce the opening hours 2. Clean the filter screen. 3. Check the temperature of the laboratory
	The temperature of the incubator is too low	
Detection stops for 30 minutes		If System is open or the power failure lasts for over 30 minutes, all the culture bottles should pass the screening.
Data loss		1. Copy the data and address the distributor. 2. All the bottles should pass the screening.

Culture bottle error	1. The bottle is improperly installed in the position 2. The bottle is unloaded from the device 3. Wrong position No. 4. Position and control device is broken.	1. Place the culture bottle properly 2. Get back the missing culture bottles to the gaging position. 3. Take out the bottle from the error position and restart System. 4. Turn off the damaged bottle position.
System cannot be energized	Electric fuse blowing	Change the 5A 230V fuse
System is full		Take out the bottles with the expired incubation time.
The rotating speed of the rotor is overrun	The rotating speed of the rotor is slow down	1. Confirm the placement of the culture bottle is correct 2. Remove the reason of hindrance 3. All the culture bottle are rejected.
	The rotor stops running	

 If it is not possible to correct the fault yourself, address the service center of the manufacturer's authorized official in the Russian Federation

5 Transportation and storage

5.1 Transportation

The packed medical devices are transported by all types of the enclosed transport according to the haulage rules applied for this type of transport.

Placing and fixing the boxes with the medical devices to provide their stable position, excluding the possibility of the boxes displacement and their clashes and impacts with the walls of the transport facilities.

Transportation conditions:

Environment temperature rate: from minus 20°C to plus 55°C;

Relative humidity range: 45% - 95%;

Atmospheric pressure range: 50.0 kPa–106.0 kPa.

5.2 Storage

When packing into the wooden boxes the allowable height of decking should be indicated on the packaging places, but no more than 2 (two).

The medical devices should be installed and stored in the following conditions to provide the normal performance:

- Environment temperature: from plus 10°C to plus 30°C;
- Relative humidity: 45- 85%;
- Atmospheric pressure: 86.0 kPa -106.0 kPa
 - The device should be kept away from intrusion of dust, direct sunlight, liquids, corrosive gases, vibration and high electromagnetic disturbance.

6 Warranty policy

The manufacturer guarantees the medical device compliance with the regulatory requirements and manufacturers' documentation in case of observance of terms of transportation, storage and operation.

Operating life – 10 years.

Guaranteed service life – 1 year since the purchase.

Guaranteed storage life – 2 years.

The manufacturer is not responsible for any inflicted damage or injuries when operating the System until his fault is proven.

7 The medical device utilization and elimination procedure

The product is subject to disposal in accordance with the current legislation of the Russian Federation at the time of disposal.

The components of the System belong to class A "Epidemiologically safe waste, the composition is close to solid household waste" and are disposed of in accordance with the requirements of SanPiN 2.1.7.2790-10.

Biological samples may contain pathogenic microorganisms, including hepatitis b virus and human immunodeficiency virus (HIV), so the media in the vials belong to class B - "Epidemiologically hazardous waste" and require appropriate precautions for their disposal according to SanPiN 2.1.7.2790-10

Used bottles are subject to disinfection. If the laboratory is equipped with an autoclave, the disinfection is performed by autoclaving (exposure to water vapor under pressure) at a temperature of 121°C and a pressure of 1.2 ATM. 20 minutes. Bottles with the remains of blood or other biological samples, closed with lids, are

collected in containers or plastic bags, suitable for the physical method of disinfection, resistant to liquid leakage, and placed in an autoclave. After autoclaving, the waste is moved to yellow bags labeled "class B Medical waste". The person responsible for the collection of waste in this medical unit ties the bag or closes it with tags or other devices and indicates on the label the date, time, Name Of the responsible person. Packages are removed to temporary storage of waste for subsequent removal and disposal. If the laboratory is not equipped with an autoclave, disinfection is carried out using other methods established by local rules and regulations.

8 Appendix A

The information about the examination with the use of American typical culture collection (ATCC) strains.

Strains	Strain code	Conditions		Pathogenic microorganism detection method	Results
<i>Neisseria meningitidis</i>	ATCC13090	Temperature (C)	Time (h)	The culture detected by System	Triggering the alerting system
<i>Haemophilus influenzae</i>	ATCC19418	35	<72	The culture detected by System	Triggering the alerting system
<i>Streptococcus pneumoniae</i>	ATCC6305	35	<72	The culture detected by System	Triggering the alerting system
<i>Streptococcus pyogenes</i>	ATCC19615	35	<72	The culture detected by System	Triggering the alerting system
<i>Pseudomonas aeruginosa</i>	ATCC27853	35	<72	The culture detected by System	Triggering the alerting system
<i>Candida albicans</i>	ATCC18804	35	<72	The culture detected by System	Triggering the alerting system
<i>Staphylococcus aureus</i>	ATCC25923	35	<72	The culture detected by System	Triggering the alerting system
<i>Escherichia coli</i>	ATCC25922	35	<72	The culture detected by System	Triggering the alerting system
<i>Alcaligenes faecalis</i>	ATCC8750	35	<72	The culture detected by System	Triggering the alerting system
<i>Clostridium histolyticum</i>	ATCC19401	35	<72	The culture detected by System	Triggering the alerting system
<i>Clostridium perfringens</i>	ATCC13124	35	<72	The culture detected by System	Triggering the alerting system
<i>Bacteroides fragilis</i>	ATCC25285	35	<72	The culture detected by System	Triggering the alerting system
<i>Bacteroides vulgatus</i>	ATCC8482	35	<72	The culture detected by System	Triggering the alerting system

9 Appendix B

Safety requirement, test method and inspection rule (according to the GB4793.1-2007 standard)

B0 Application condition of the instrument

- a) For indoor use.
- b) Temperature: 10-30°C
- c) Humidity: No more than 85%
- d) Power condition: The power fluctuation is no more than 10% of the standard voltage
- e) Overvoltage category: The transient overvoltage category II
- f) Pollution level: Grade II

B1 Mark and file

B1.1 Mark

B1.1.1 Overview

Requirement: The words and symbols marked outside the device should be conformed to the requirement of 5.1.1 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.1.2 Mark

Requirement: The mark indicated on the device should meet the requirement of 5.1.2 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.1.3 Mains supply

Requirement: The mark indicated on the device should meet the requirement of 5.1.3 in GB4793.1-2007.

- a) Nature of supply.
- b) The rated value of the supply voltage.
- c) The maximum rated power.

Test method: Check the quality using visual inspection.

B1.1.4 Fuse

Requirement: It should meet the requirement of 5.1.4 in GB4793.1-2007.

Test method: Check the quality using visual inspection.

B1.1.5 Terminal, connector and operation equipment

Requirement: The purpose of the terminals, connectors and indicators should be specified if it is necessary for safety. If there is no enough space, the symbol 14 of Table 1 can be used.

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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ShaoYunjie

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(Date: Sep. 23, 2019)

证明书查询 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

APPROVED BY

SCENKER Biological Technology
Co., Ltd.

CUI GE/Legal

representative

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September 12, 2019

Automated bacteriological system ЮHOHA® LABSTAR 100
with the accessories

SERVICE MANUAL

Preface

We thank you for buying the Automated bacteriological system LABSTAR 100.

LABSTAR 100 provides the rapid growth of microorganisms and is designated for the clinical culture and detection of the pathogens in the blood and other aseptic body fluids. The system allows quickly and accurately detecting the pathogenic microorganisms in the blood of the patients with septicemia, bacteremia, etc. Besides, this device can also be used to detect the pathogenic microorganisms in the body cavities that are usually without bacteria (such as thoracic cavity, abdominal cavity, articular cavity, pericardial cavity and cerebrospinal cavity, etc.), that allows making the clinical diagnosis quickly and prescribing the effective treatment regimen.

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The System is recommended to be used in the laboratory that meets the requirements of Level-2 bio-safety laboratory. During the device operation, it is necessary to take the bio-safety controls. The medical waste, such as the discarded culture bottles, must be sterilized with the high-pressure steam or other harmless treatment method according to the local regulations.

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The instrument will be presented under **ЮЮЮЮ**[®] trademark on the Russian market.

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Revision date: December 01, 2014

The manual revision for the registration in Russian Federation:

v.2 dated September 12, 2019.

Scenker Biological Technology Co., Ltd.

Contents

1 Medical device data	5
1.1 Medical device name	5
1.2 Information about the Medical device manufacturer	5
1.3 Information about the Medical device manufacturer authorized official	5
2 Operation and description	5
2.1 Medical device overview	5
2.2 Technical details, parameters and characteristics of the medical device	6
2.3 The contents (contents of delivery) of the medical device	23
2.4 Medical device precaution measures and limits of use	23
2.5 The description of the package, marking and warning signs on the device	25
3 Medical device operation	28
3.1 Environmental requirements	28
3.2 Installation procedure	29
3.3 Setting of system parameters	32
3.4 Operation procedure	37
3.5 Culture media and application scope	45
3.6 How to collect the blood and other sterile biological fluid samples	45
3.7 Results explanation	48
4 Maintenance and repairs	49
4.1 The user-performed technical maintenance:	49
4.2 The common failures and treatment measures	50
5 Transportation and storage	51
5.1 Transportation	51
5.2 Storage	52
6 Warranty policy	52
7 The medical device utilization and elimination procedure	52
8 Appendix A	54
9 Appendix B	55
10 Appendix C	66
11 Appendix D	68
12 Appendix E	71

1 Medical device data

1.1 Medical device name

Automated bacteriological system **ЮHOHA®** LABSTAR 100 with the accessories (hereinafter referred to as System).

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Concerning the quality and service of the Medical device “Automated bacteriological system **ЮHOHA®** LABSTAR, design versions: **ЮHOHA®** LABSTAR 50, **ЮHOHA®** LABSTAR 100 with the accessories”, manufactured by SCENKER Biological Technology Co., Ltd., China, please contact the manufacturer’s authorized official in Russian Federation:

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- Culture Medium Bottles **ЮHOHA®** for aerobe cultivation (color of label and removable cover is red, designation ■) – a standard culture medium bottles designed for adult patients not treated with antibiotics, the volume of the biological sample is from 5 ml to 10 ml;

- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for aerobe cultivation (color of label and a removable cover is blue, designation ■), culture medium bottles designed for adult patients treated with antibiotics, the volume of the biological sample is from 5 ml to 10 ml;

- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for anaerobe cultivation (color of label and removable cover is cadmium orange, designation ■), culture medium bottles designed for adult patients treated with antibiotics, the volume of the biological sample is from 2 ml to 10 ml;

- Antibiotic-neutralized Culture Medium Bottles **ЮHOHA®** for aerobe cultivation for children (the color of the label and removable cover is pink, designation ■), culture medium bottles designed for children, the volume of the biological sample is from 2 ml to 10 ml.

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
System ЮHOHA® LABSTAR 100 consists of (fig.1):

- Main module;
- Main module power cord;
- Computer with the preinstalled software;
- Data cable;
- Computer power cord;
- Barcode scanner (if necessary).

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*Fig. 1 – Automated bacteriological system ЮHOHA® LABSTAR 100 with the accessories
1 – Main module, 2 – Computer, 3 – Monitor, 4 – Keyboard, 5 – Mouse, 6 – Barcode scanner,
7 – Power indicator, 8 – Temperature indicator, 9 – Software*

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LABSTAR 100 is designed as a bracket consisting of five lanes of rotating frames, each lane includes 20 bottle Position Nos. The bracket size reaches $(475 \times 545 \times 76) \pm 1$ mm, the diameter of the bottle Position Nos. is $36.5 (+0.1-0.3)$ mm. Each Position No. is connected to a highly sensitive fluorescent sensor. The total capacity is 100 bottles. The Position Nos. are numbered from 1 to 100: 001-100. The operation cycle (scan mode) of System occurs every 10 minutes..


Fig.2 - LABSTAR 100 System turning lane bracket

The analytical unit of the system computer creates the growth curve based on the measurement data, and its form allows making the conclusion of the presence or the absence of the pathogenic microorganisms in this sample.


Fig.3 –LABSTAR 100 main screen. Growth curve

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
Fig. 8 – Main control PWB


Fig. 9 - LABSTAR 100 System process functional diagram

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Power adapter capacity	AC220-240V $\geq$ 200W
Type/Quantity	USB 2.0 type A/4, USB 3.0 type All-WO/1, PS/2/2, RS232, RJ45, VGA
Disc drive type	DVD
Software	Preinstalled, LABSTAR 100, version 3.0, Copyright certificate of registration No. RZDZ 0543399
Overall dimension, mm	155x400x370

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Software characteristics	Model LABSTAR 100
Name	LABSTAR 100
Version	V3.0
Platform	Intel ^R Celeron ^R G1840
Operating system	Microsoft Windows Professional10
Software programming language	Delphi
Software data transmission communication interface	RS-232
Software data transmission communication interface	Automatic
Software data storage communication interface	-
Software response at full load, s, max	1
Accident recovery time (complete software failure due to a technical or operational accident), h, max	2
Processing speed, Mbit /s	100
Recovery time, if virus programs are detected, h, max	2
Performance, Gb/s, min	1.5

Search time at full load, s, max	1
Interfaced systems	Ethernet(TCP/IP)
Protection means	Software account
Safety class	P "Minor injuries are possible"


Fig.11 - LABSTAR 100 screen with the system information

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
Fig. 12 – LABSTAR 100 Main computer interface
1-Bottles position, 2- control panel, 3- growth curve,
4- system information

(1) – The position of the bottles on the turntable is always displayed, the empty positions have the white color, the loaded positions are yellow, the positive bottle

positions are red and the negative bottle positions are green. The error positions are blue.

(2) - The control panel located in the screen footer consists of buttons:

Bottle load	To load the bottles
Bottle removal	To remove the bottles
Sample data	To enter data on the patient
Growth curve	To view a graphic representation of the curve
Statistics	To request data
System settings	To set up the system
About the system	To receive information on the system

(3) – When passing the mouse to the loaded position, the diagram displays the growth curve. The data readout takes place each 10 minutes. The form of the curve reflects the presence or the absence of the microorganism growth in the loaded bottles. By default, the range of Y-axis is 700-4,000, the range of X-axis is 0 - the maximum time of cultivating this bottles in terms of the days. If the bottle is defined as the positive one, the chart shows the “alarming point” - the detection point and the cultivating time before the bottle positive definition. The scale of X and Y axes can be changed using the right mouse button. The date and time of the bottle installation and the date and time of alerting are pointed out under the growth curve.

(4) – System information:

- The total number – the total number of the loaded bottles;
- Under test – the number of the bottles under test;
- Positive – the number of the positive bottles;
- Negative – the number of the negative bottles.
- Temperature – culture/incubation temperature;
- State of the system: “normal” at the normal temperature, “cooling” – if the temperature is over the limit, “heating” – if the temperature is below normal;
- Warning notice "Scanning, do not open the door" - when reading the data

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Parameter	Characteristic
Model	LS-1881
Manufacturer	Guangzhou OWL Electronic Technology Co. , Ltd
Projector	Laser

Scanning frequency	150 scan/sec
Barcode contrast	$\geq 30\%$
Wave-length	650 nm
Operational temperature	-45~45°C
Resolution	PCS 0.9 (0.1~0.625mm)
Focus	10mm~420mm
Indication (LED, sound etc.)	Double color LED
Protection level	IP54
Data transmission rate	200 characters per second
The mean of communication of the scanner with System;	RS232 serial data transmission
Readable codes	Single two-dimensional code

Fig. 14 - Barcode scanner

2.2.4.7 System parameters and characteristics

Table 5. Index comparison table

Index	LABSTAR 100
Bottles capacity	100
Culture temperature	35±1.5°C
Rotating speed	18 r/min
Operation cycle (scan mode)	every 10 minutes
Rotating mechanism type	Rotating frame bracket
The system is designed for microorganisms culture and detection	+
The rotating and shaking of the bottles is provided	+
Automatic correction and quality control	+
The number of positions for the bottles	100

The availability of the graphic interface	+
The availability of the Barcode scanner	+
Fluorescent method (instrumental estimation), bottle bottom color change (manual estimation)	+
Incubation conditions: fixed temperature, rotating and vibration	+
A/C power source, V	220 ± 22
Current frequency, Hz	50 ± 1
Input power, max VA	450
Main module external dimensions, mm	$(620 \times 650 \times 800) \pm 10$
System main module weight, kg	110 ± 2
Efficiency*	20 samples a day
Time of cultivation	120 hours
Volume of biological sample	5-10 ml aerobic bottles and 2-10 ml of anaerobic and children

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2.2.4.8 Basic requirements

- a. Operating conditions
 - Environment temperature: 10-30°C;
 - General humidity: 45-85%;
 - Atmosphere pressure: 86.0-106.0 kPa;
 - Electric power conditions: 220V $\pm$ 22V, 50Hz $\pm$ 1Hz;
 - Input power: $\leq$ 450VA;
 - Away from the sources of the high electromagnetic disturbance and bright directed light;
 - System should be dust-proof; the air should not contain the acids, salts and other corrosive gases;
 - The proper grounding should be provided;

- It is necessary to control the vibration at the workplace.
- b. Appearance
 - The equipment should be flat, clean and of the uniform color, it should be without the scratches, blinks and mechanic damages.

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- The culture temperature should remain within the range of $35 \pm 1.5^{\circ}\text{C}$.
- The temperature variations should not exceed 3.0°C .

e. System functions

- Continuous culture.
- Automated monitoring and determination of the positive and negative results.

- Immediate notification about the positive bottles detection.
- Alarm signaling in case the temperature exceeds the limits of the acceptable range.

f. Positive results repeatability

The results of the testing using the standard strains of the American typical culture collection (ATCC) should be positive (see Appendix A).

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Characteristics	LABSTAR 100
Brand, model	LENOVO, model LT2024WA
Photo	
Screen size	20.0"
Matrix type	TFT TN
Screen format	16:9
Resolution	1,600x900
Brightness	250 kcd/m ²
Contrast ratio	1000:1
Response time	5 ms
Connector type	VGA, DVI
Power cable	Chinese Standard: CPCS-CCC Model: PSB-10/ST3 Type I connector, 3 flat contacts. Load-carrying capability: 10A/250B Length, min: 1.5 m

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Table 7. Keyboard technical characteristics

Characteristics	LABSTAR 100
Description	Standard classical desktop corded keyboard with 104 keys
Cable length, min, m	1.5
Connection type	USB

Photo	
-------	--

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Characteristics	LABSTAR 100
Description	Optical corded mouse
Cable length, min, m	1.5
Connection type	USB
Photo	

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- Main module – 1 pcs.
- Data cable - 1 pcs.
- Main module power cord – 1 pcs.
- The PC with the preinstalled software – 1 pcs.
- Computer power cord – 1 pcs
- Barcode scanner – 1 pcs. (if necessary)
- Service manual – 1 pcs.
- Certificate – 1 pcs.

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 During the operation of System the operator can be subjected to dangerous microorganisms, so it is recommended to install the device in laboratory premises that meet the requirements of biosafety level 2. The provisions should be made for the biological protection measures and it is necessary to fully comply with the requirements for the organization of the microbiological research laboratories. The used bottles and other medical wastes must be disposed of by autoclaving sterilization or other methods determined by the local rules and regulations (see 7 The medical device utilization and elimination procedure).

 There is the danger of infecting during testing the samples with System. As a result, the marking in the form of the Biohazard Mark is provided on the device.

 It is prohibited to touch the electric power supply of the operating System.

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 The device back panel has the mark “Caution! Please, refer to the service manual”.

 In case of the noncompliance with the operating procedure established in this document, the normal operation of System can be disrupted. In this case the findings can be distorted, and the destruction of the equipment as well. The user and (or users) of the device hold responsibility of the consequences of this noncompliance.

 System main module has the drive unit connected to the A/C power source with the voltage of 220 V and it can be a source of the external electromagnetic disturbance. (see Appendix D EMS)

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- Insufficient sample volume: within the scope of the sample volume (Adults: 5-10 ml, children: 2-10 ml) proportionally to the positive growth of the cultivated microorganisms.

- The sample collection within 24 hours is not enough: the bacteremia is remitted in 85% cases, several samples collection can assure the positive growth recognition.

- The blood sampling should be done within one locus every time.
- The inappropriate time of the blood sampling.
- Carrying out the antimicrobial therapy before or during the blood sampling.

- Storing the samples at the room temperature : it is not allowed storing the samples in the freezer, because the bacteria with the special nutritive demands die at the low temperatures.

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- Power supply conditions, input power;
- Product standard code;
- Production number;
- Production date;
- Manufacturer's name, address, telephone;
- Authorized official's name, address, telephone;
- Product license number and date,

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

- Manufacturer's name and address;
- Distributor's name and address;
- Product name, model;
- Serial number;
- Weigh net/ gross;
- Volume (length x width x height);
- Manufacture date;
- Product standard code;

- Indication “Fragile, handle with care” (the corresponding symbol is marked);
- Indication “Keep dry” (the corresponding symbol is marked);
- Indication “Up” (the corresponding symbol is marked);
- Indication “Do not roll” (the corresponding symbol is marked);
- Indication “Limit of the pile levels number”;
- Indication “Fragile precision tool. Caution”;
- transportation and storage conditions (the corresponding symbols are marked);

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3.1 Environmental requirements

3.1.1 Operational conditions

To provide the normal operation System should be installed and stored in the following environment:

- Environmental temperature: 10°C-30°C;

- Relative humidity: 45- 85%;
- Atmospheric pressure: 86.0 kPa ~106.0 kPa
- A/C power source: 220±22 V, frequency 50±1 Hz;
- Keep System away from intrusion of dust, bright directed light, liquids, corrosive gases, vibration and high electromagnetic disturbance

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3.1.3 Power requirements

The analyzer is designed for the AC power supply with the voltage of 220 V ± 22 V, the frequency of 50 Hz ± 1 Hz; The device input power is 450 VA. There must be the regulated power supply between the device and the power system. The electric fuse characteristics are indicated on the back panel. When replacing the electric fuse, it is important to take the power cable characteristics into account.

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3.2.1 Unpacking and acceptance

Open the package, carefully remove System with the accessories from the packing box. System packing list is inside of the box. After the unpacking, please check the contents of delivery according to the packing list. In case of the

deficiency, address the distributor immediately.

3.2.2 Installation and implementation

The preoperational phase should be carried out in accordance with the requirements listed in this section.

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
Fig.19 - Back panel of LABSTAR 100 main module
 1-4 – Fans, 5 – Switch, 6 - power connector + electric fuse,
 7 - RS232 port for connecting the data cable

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The registration process that is carried out by the manufacturer's authorized official is required for the initial system start up. In order to do that, enter the organization name in "Department name" field and the System code indicated in the delivery set accompanying documentation, into "Device code" field of the "Logging in" pop-up window. Then obtain the license number from the manufacturer and enter it in the "License" field.


Fig.21 - Logging in

After pressing the "Confirm" button the device will be registered and ready for use.

3.3 Setting of system parameters

3.3.1 Culture time

The system automatically sets the culture time for 120 hours. If there are special requirements, the time parameter on the screen should be modified when the bottles are put in and the time set on the "Set time" screen (p. 3.3.3.2).


Fig. 22 – Default culture time is 120 hours

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3.3.3 System settings

3.3.3.1 Department name

Enter the instrument control procedure, click “System settings” and the interface of the dialog box appears. Click “Department” and then click “Add” to input the name of the specific department. Click “Delete” to delete the chosen department.


Fig.23 - “System Settings” menu, department name

3.3.3.2 Signal mode

If necessary, an acoustic signal informing about positive or negative bottle can be adjusted.

Enter the instrument control procedure, click “System settings” and the interface of the dialog box appears. Click "Signal Mode", then select an audio file from a designated folder and confirm the selection.


Fig.24 - "System Settings" menu, signal mode

In this window, the culture / testing time can be set, by default it is 120 hours (section 3.3.1).

3.3.3.3 Bottle

In this window, the system settings can be reset by clicking "Remove all samples" button. This will lead to the bottle position and testing results data reset.

⚠ Use the option of the system settings reset only in exceptional cases, for example, for serious instrument malfunctions or power cutoff for more than 1 hour. It is not recommended to use this option for normal operation of the instrument, otherwise all the data will be lost!!!

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
Fig.25 - "System Settings" menu, bottle

This window also allows marking the defective Position No. to avoid loading bottle in it.

This requires entering the number of the defective Position No. (for example, Position No. 2) in the right field and click "Confirm" button. This Position No. will be displayed in black.


Fig.26 - LABSTAR 100 screen with a defective Position No.

After troubleshooting, click "Delete selected parameters". "o" will appear in the Position No. field and the Position No. will be white, i.e. empty (Fig. 27).


Fig.27 - Recovery of the Position No. parameters

3.3.3.4 Database

The window allows for configuration:

- Deleting data from the test log;
- Manual backup of data;
- Connection to the LIS system.


Fig.28 - "System Settings" menu, database

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Blood and other body fluids sample collection should be carried out according to the rules adopted in the laboratory, using the special bottles manufactured by SCENKER Biological Technology Co., Ltd. (see the more detailed information about the analysis of blood and other body fluids in the Cultural medium Bottles documentation).

- The blood and other body fluids sample for the culture should be taken with keeping the rules of the aseptic to avoid the microorganisms' penetration from the outside environment.
- Before using the bottles, medium transparency should be visually assessed, any opacity means their unfitness.
- The bottles should be warmed for the room temperature (18-30°C).
- The culture should be carried out during the temperature rise, when the fever begins. To take blood using special sterile systems for blood sampling or sterile syringes in the bottle created a vacuum, which facilitates the sowing of blood samples in the bottle. Break off the plastic covers from the bottles and sterilize the inside cover with the drape saturated in 70% ethanol (or other substitutive antiseptic) within one minute. Carefully sterilize the skin over the needled vein with 70% ethanol and then with 1-2% tincture of iodine within 30 seconds. After drying of the treated area, perform venipuncture. For adult patients it is recommended to use a pair of aerobic and anaerobic bottle. When using the system for blood sampling, the inoculation is performed first in an aerobic vial and then in an anaerobic one. (Do not remove the needle from the patient's vein during blood collection. After the end of the first extract needle from vial, then from veins

patient). When using a syringe, first inject blood into an anaerobic vial and then into an aerobic one (Avoid air entering the anaerobic bottle). The volume of the blood is from 5 ml to 10 ml for the adults for aerobic bottle and from 2 ml to 10 ml for the anaerobic and children bottle (the marking on the name plate of the bottle allows you to estimate the volume of the taken blood). Optimal pathogen isolation is achieved by adding the maximum sample volume. Mix the bottle contents carefully with the circular motions. After finishing the procedure, the iodine should be removed from the skin with the drape using the 70% ethyl and the puncture area should be covered the adhesive plaster. The bottles should be marked with indicating the name, initials, date of birth and sampling date, and then transported to the laboratory in the at the room temperature (18-30°C) in a dark place (not in the freezer) for no more than 24 hours.

- Cerebrospinal fluid is obtained with lumbar puncture, less often with puncture of the lateral ventricles of the brain, preferably before taking antibiotics. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. The first portions of liquor have the highest potential for contamination by microorganisms from the skin and should be discarded. Immediately after the puncture, not less than 5 ml and not more than 10 ml of the resulting liquor should be aseptically introduced into an aerobic vial (in certain clinical situations, an anaerobic vial is required). Large volumes (5-10 ml) increase the sensitivity of the culture and are necessary for optimal isolation of whimsical bacteria and fungi. If the volume of liquor is small from 2 to 5 ml, as well as liquor from children - you should use aerobic bottles for children. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of pathogens. During storage and transportation, excessive overheating and especially cooling should be avoided (a thermos or a thermal container with a temperature of $35 \pm 1.5^{\circ}\text{C}$ can be used). Cooling the liquor below 30°C leads to the death of meningococci. Storage and transportation time should not be more than 8 hours.

- In the study of body fluids from normally sterile cavities (peritoneal, synovial, pleural, articular, pericardial fluid) for adults, it is recommended to take for each aerobic vial from 5 ml to 10 ml, for each anaerobic vial from 2 ml to 10 ml of the test liquid, for children for each vial it is necessary from 2 ml to 10 ml of the test liquid. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. Immediately after the puncture, the biological fluid should be aseptically introduced into the aerobic vial (in certain clinical

situations, an anaerobic vial is required). The minimum volume of liquid for bacteria isolation is 2-5 ml, for fungi isolation-10 ml. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of microorganisms. Store at room temperature (18-30°C) for up to 24 hours.

- The rate of growth of pathogens is directly proportional to the volume of biomaterial. Optimal pathogen detection is achieved with maximum sample volume.
- After sampling and inoculation of blood or biological fluids, it is necessary to tear off half of the label with a bar code from the bottle (the other half of the label can not be torn off) and stick it on the form of the referral for laboratory research


⚠ When separating the label with the barcode, it is necessary to check that the information on the issue is the same as on the bottle!

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3.4.2 Bottles loading

3.4.2.1 Using the barcode scanner

Open the door of the main module. Click “Load”. The dialog box will automatically pop up on the screen. Use the barcode scanner to read the barcode from the label. System will automatically allocate the Position No. for placing the bottle. Put the bottle to the necessary position and click “Load” (at the same time the Position No. at main program interface manifests yellow) (fig.29). Load all the bottles for the check and click “Confirm”. Close the maid module door and it will automatically resume operation.


Fig. 29 - Loading bottles into Position No. 66

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Information Input: Enter the control procedure of the instrument, click “Data”, and the interface of dialog box will appear.

Fig. 30 - Sample data entry and request

Enter the sample number in the appropriate field using a special barcode scanner or manually. Following that, fill the patients data (Surname, first name & patronymic, gender, age, department name, patient number, patient history No., date of check, sample category). Click "Save" after patient data is filled.

To request data, it is necessary to enter the sample number and click "Search", the system will provide all information about the sample.

To receive information in a graphical form, enter the instrument control procedure and click "Growth Curve" and the interface of dialog box will appear. Click "Request" and the interface of request dialog box will appear. Use a dedicated barcode scanner. The instrument will automatically display patient data, study parameters and results, and the growth curve.


Fig. 31 - Patient data request, growth curve

Statistics: Enter the instrument control procedure, click "Statistics" and the interface of statistics dialog box will appear. Statistics can be filtered by date, department, patient's name, study results, sample number.


Fig. 32 – Statistics by patient's name


Fig. 33 - Statistics by date

3.4.4 Culture

After the culture bottle is put into the main module, System will automatically start the culture. The temperature should be maintained at the level of $35 \pm 1.5^\circ\text{C}$. The LABSTAR 100 turntable will rotate with the constant speed of 18 r/min.

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3.4.6 Alerting system

The system will analyze the detection signal. It will alarm immediately when it discovers the positive bottle. The Position No. of the culture medium in the main program interface manifests red after the alarm. Unload the bottle from the device and use the sample for the subculture to conduct the additional analysis experiments. The culture bottle with no growth of culture medium of microorganism after 5 days of successive culture will be reported as the negative bottle and the corresponding alarm signal will occur. The Position No. of the culture bottle in the main program interface will manifest grass-green and the culture bottle can be taken out.


Fig. 34 - Main program interface with red and green bottles

After all the positive bottles are taken out, the sub-culture should be applied immediately and the additional analysis should be carried out so as to prevent the instrument from alarming falsely!

3.4.7 Unloading the bottles

Open the device door and unload the negative bottles (green color on the interface screen) or the positive bottles (red color on the interface screen). Click "Unload", the instrument will automatically pop up the dialog box, then use the bar code scanner to read the sample bar code. In the meanwhile, the Position No. of the culture bottle in the main program interface manifests white, which indicates that the culture bottle has been taken out. Unloading bottles by sample number is possible using the barcode scanner or entering the number manually. Remove all bottles ready for unloading and click "Confirm".

Fig. 35 - Unloading the bottle

After the culture bottle is taken out, please immediately close the instrument door, so as to prevent low temperature inside the instrument!

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
Fig. 36 LABSTAR System operation cycle flow diagram

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3.5.1 Culture media

Table 9

Name	The volume of the collected material	The volume of the medium	Application	Barcode/ ID No.	Photo
Culture medium Bottles ЮHOHA[®] for aerobic cultivation:	5-10 ml	25 ml	The samples of the patients who did not take the antibiotics	XBZ	
Antibiotic-neutralized Culture medium Bottles ЮHOHA[®] for aerobic cultivation:	5-10 ml	25 ml	The samples of the patients who took the antibiotics	XCR	
Antibiotic-neutralized Culture medium Bottles ЮHOHA[®] for anaerobe cultivation	2-10 ml	25 ml	The samples of the patients who took the antibiotics	XYX	
Antibiotic-neutralized Culture medium Bottles ЮHOHA[®] for aerobic cultivation for children.	2-10 ml	25 ml	The samples of the children or other patients that are difficult to collect	XET	

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3.6.1.4 Skin sterilization process

- Sterilize the skin over the needled vein with 70% alcohol solution within 30 seconds;
- Sterilize with tincture of iodine within 30 seconds or tamed iodine within 60 seconds;
- Remove the rests of iodine from the skin with 70% alcohol solution.
- Carry out the blood collecting after the alcohol evaporation.

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➤ If it is impossible to send the sample for the examination at once, it can be stored at the room temperature or in the constant-temperature bath (35-37°C). Cold storage is not allowed;

➤ The delay for the examination cannot exceed 20 h after collection.

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➤ The bottles are rejected if they are kept for more than 20 h after the collection, the bottle type does not correspond to the examination type, the bottle is expired, the volume is not enough for the examination etc.

3.6.2 Cerebrospinal fluid is obtained with lumbar puncture, less often with puncture of the lateral ventricles of the brain, preferably before taking antibiotics. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. The first portions of liquor have the highest potential for contamination by microorganisms from the skin and should be discarded. Immediately after the puncture, not less than 5 ml and not more than 10 ml of the resulting liquor should be aseptically introduced into an aerobic vial (in certain clinical situations, an anaerobic vial is required). Large volumes (5-10 ml) increase the sensitivity of the culture and are necessary for optimal isolation of whimsical bacteria and fungi. If the volume of liquor is small from 2 to 5 ml, as well as liquor from children - you should use aerobic bottles for children. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of pathogens. During storage and transportation, excessive overheating and especially cooling should be avoided (a thermos or a thermal container with a temperature of $35 \pm 1.5^{\circ}\text{C}$ can be used). Cooling the liquor below 30°C leads to the death of meningococci. Storage and transportation time should not be more than 8 hours.

3.6.3 In the study of body fluids from sterile normal cavities (peritoneal, pleural, articular, pericardial fluid) for adults is recommended to take for each aerobic vials from 5 ml to 10 ml, for each anaerobic vial of 2 ml to 10 ml of the test liquid, for children of each vial must be from 2 ml to 10 ml. The puncture site is treated with 70% ethyl alcohol, then with 1-2% tincture of iodine within 30 seconds. Immediately after the puncture, the biological fluid should be aseptically introduced into the aerobic vial (in certain clinical situations, an anaerobic vial is required). The minimum volume of liquid for bacteria isolation is 2-5 ml, for fungi isolation-10 ml. Samples should be sent to the laboratory for examination as soon as possible. This will increase the probability of detection of microorganisms. Store at room temperature ($18-30^{\circ}\text{C}$) for up to 24 hours.

3.7 Results explanation

The human blood usually has no bacteria. The bacterial invasion can cause the bacteremia or septicemia. The patients with endocarditis, endarteritis, early typhoid and brucellosis are considered to be suffering from the bacteremia. Some other infections may be accompanied with intermittent bacteremia. Various results of the blood culture should be analyzed according to the specific situation. The culture results do not provide the exact bacteremia indication.

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- (1) There is no bacteremia
- (2) The patient has the clinical signs but the bacteremia was not detected

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Type of the maintenance	LABSTAR 100
Cleaning the air filters with the fresh water	Once per 2 weeks
Dust collection in System functional area	Once per 2 months
Cleaning the liquids spilled on System	In case of the leakage or the bottle damage, remove the spilled liquid according to the decontamination procedure recommended in the laboratory. After the disinfection, clean with the wet wipe (moistened with water only) and rub dry.

Cleaning the liquids spilled inside of System	Remove the samples spilled inside of System and carry out the disinfection according to the authorized procedure.
Checking System cleanness	Once per 3 months, clean with the dry cotton plug if necessary
Checking the output voltage of the regulated power supply	Once per 6 months
Fixed temperature and humidity maintenance in the laboratory	Within 24 hours
System supply disconnection in case of the power outage	Turn on System after the power restart
Switch off System power supply and restart it in case of the sputter of the signal system and failure to correct the error	In 3 minutes after shutdown

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Error category	Cause	Solutions
Main module power fail	No power cable or the power cable unplugged.	Check the power cable connection at both ends.
	The main module power switch is in the “off” position.	Set the power switch into “On” position (it is located in the lower right corner on the main module back panel).
Lost connection between the PC and the main module	No data cable between the PC and the main module or it is wrongly connected.	Check the cable connection at both ends. Restart the main module.
Temperature Alarm	The temperature of the incubator is too high	1. Reduce the opening hours 2. Clean the filter screen. 3. Check the temperature of the laboratory
	The temperature of the incubator is too low	

Detection stops for 30 minutes		If System is open or the power failure lasts for over 30 minutes, all the culture bottles should pass the screening.
Data loss		1. Copy the data and address the distributor. 2. All the bottles should pass the screening.
Culture bottle error	1. The bottle is improperly installed in the position 2. The bottle is unloaded from the device 3. Wrong position No. 4. Position and control device is broken.	1. Place the culture bottle properly 2. Get back the missing culture bottles to the gaging position. 3. Take out the bottle from the error position and restart System. 4. Turn off the damaged bottle position.
System cannot be energized	Electric fuse blowing	Change the 5A 220V fuse
System is full		Take out the bottles with the expired incubation time.

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Used bottles are subject to disinfection. If the laboratory is equipped with an autoclave, the disinfection is performed by autoclaving (exposure to water vapor

under pressure) at a temperature of 121°C and a pressure of 1.2 ATM. 20 minutes. Bottles with the remains of blood or other biological samples, closed with lids, are collected in containers or plastic bags, suitable for the physical method of disinfection, resistant to liquid leakage, and placed in an autoclave. After autoclaving, the waste is moved to yellow bags labeled "class B Medical waste". The person responsible for the collection of waste in this medical unit ties the bag or closes it with tags or other devices and indicates on the label the date, time, Name Of the responsible person. Packages are removed to temporary storage of waste for subsequent removal and disposal. If the laboratory is not equipped with an autoclave, disinfection is carried out using other methods established by local rules and regulations.

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Not applicable.

11 Appendix D

Automated bacteriological system **ЮHOHa®** LABSTAR, design versions:

ЮHOHa® LABSTAR 50, **ЮHOHa®** LABSTAR 100 with accessories meets the requirements for emissions and noise immunity GB/T 18268.1-2010 and GB/T 18268.26-2010.

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Automated bacteriological system **ЮHOHa®** LABSTAR, design versions:

ЮHOHa® LABSTAR 50, **ЮHOHa®** LABSTAR 100 with accessories designed for use in basic electromagnetic environment as below.

Electromagnetic emission

Electromagnetic emission test	Conformity	Electromagnetic environment - instructions
Radio emission	Group 1	Automated bacteriological system ЮHOHa® LABSTAR, design versions: ЮHOHa® LABSTAR 50, ЮHOHa® LABSTAR 100, with accessories uses radio emission only to perform its internal function. Thus, the radio emission of the product is extremely small and probably cannot cause interference electronic equipment located nearby.
Radio emission	Class A	Automated bacteriological system ЮHOHa® LABSTAR, design versions: ЮHOHa® LABSTAR 50,

		ЮHOHa* LABSTAR 100, with accessories suitable for use in all locations except residential buildings and buildings directly connected to the distribution network supplying residential buildings. When used in a residential area, it may cause radio interference. In such cases, the user should take appropriate measures.
Harmonic currents	Class A	
Voltage fluctuations and flicker	Corresponds	

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The medical device “Automated bacteriological system **ЮHOHA®** LABSTAR, design versions: **ЮHOHA®** LABSTAR 50, **ЮHOHA®** LABSTAR 100 with the accessories” corresponds to the national standards of the Russian Federation for the products:

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