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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01702936

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SIEMENS SHENZHEN MAGNETIC RESONANCE LTD. on the annexed DOCUMENT is genuine.

China Council for the Promotion
of International Trade

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验证网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

**Submission document
for Product Registration**

Legal Manufacturer: Siemens Shenzhen Magnetic Resonance Ltd.

Product: MAGNETOM Free.Max

The following document has been printed by the manufacturer in China.

- User manual supplement for Russia

Sincerely yours,

Siemens Shenzhen Magnetic Resonance Ltd.

Signature

Name

Yu Ruo Xi
(Regulatory Affairs Specialist)

Date

July 27, 2023

Manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.
Siemens MRI Center, Gaoxin C. Ave., 2nd Hi-Tech Industrial Park, 518057
Shenzhen, PEOPLE'S REPUBLIC OF CHINA

User manual supplement for Russia

*The magnetic-resonance system MAGNETOM Free.Max
with accessories*

1 The name of the Medical Device

The magnetic - resonance system MAGNETOM Free.Max with accessories.

2 Manufacturing information of the Medical Device

Manufacturer

Name:	Siemens Shenzhen Magnetic Resonance Ltd.
Address:	Siemens MRI Center, Gaoxin C. Ave., 2nd Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
Phone:	+86-755-26525421
Web:	www.siemens-healthineers.com

Manufacturing site

1) Name:	Siemens Shenzhen Magnetic Resonance Ltd.
Address:	Siemens MRI Center, Gaoxin C. Ave., 2nd Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA
Phone:	+86-755-26525421
Web:	www.siemens-healthineers.com
2) Name:	Siemens Healthcare GmbH, Magnetic Resonance (MR)
Address:	Henkestr. 127, 91052 Erlangen, Germany
Phone:	+49 9131 84-0
Web:	www.siemens-healthineers.com

Local representative in Russia

Name:	Siemens Healthcare LLC
Address:	Dubininskaja Str. 96, 115093 Moscow, Russia
Tel.:	+7 495 737 12 52
Web:	www.siemens-healthineers.com/ru/

3 Designation of the medical device

Intended use

The MAGNETOM MR system is indicated for use as a magnetic resonance diagnostic device (MRDD) that produces transverse, sagittal, coronal, and oblique cross sectional images that display the internal structure and/or function of the head, body, or extremities. Other physical parameters derived from the images may also be produced. Depending on the region of interest, contrast agents may be used.

These images and the physical parameters derived from the images when interpreted by a trained physician yield information that may assist in diagnosis.

The MAGNETOM system is not a device with measuring function as defined in the EU regulation 2017/745 (MDR). Quantitative measured values obtained are for informational purposes and cannot be used as the only basis for diagnosis.

Your MR system is a medical device for human use only!

Indications

Magnetic Resonance Imaging (MRI) is a noninvasive technique used for diagnostic imaging. MRI with its soft tissue contrast capability enables the trained user to differentiate between various soft tissues, for example, fat, water, and muscle, but can also discriminate between the different bone structures. This allows for visualization of anatomic details and can offer more sensitivity and specificity in the detection of certain abnormalities than other imaging modalities, when read and interpreted by trained personnel. In general, it is the health care provider who decides on the appropriate diagnostic procedure based on, for example, the patient's physical examination and medical history.

Intended users

The MR system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with countryspecific regulations, for example, physicians, trained radiological technicians or technologists, having received the necessary user training.

This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines and emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual.

Please contact Siemens Service for more information on available training options and suggested duration and frequency of such training.

Intended target population

As MR imaging is a general-purpose imaging modality, the system is intended for use for all humans that do not fall under the contraindications mentioned in the Contraindications chapter.

Contraindications and possible side effects

In principle, a qualified physician must evaluate the risk/benefit ratio of the MR examination for every patient.

To date, there is no scientific proof that MR examinations are harmless for pregnant women, the unborn (embryos or fetuses) and children under two years of age. An embryo or a fetus is especially sensitive to potential thermal events during the first three months of pregnancy. Therefore, avoid scanning pregnant patients especially in the first three months of pregnancy and avoid scanning patients with unknown pregnancy status.

Typical contraindications for MR examinations are:

- Electronic implants, for example, dedicated pacemakers, stimulators, insulin pumps
- Electrically conductive implants and prostheses
- Artificial heart valves, aneurysm clips
- Artificial anus (anus praeter) with magnetic closure
- Transdermal drug patches with metallic backings
- Metallic spirals for contraception (IUDs = intrauterine devices)
- Metal splinter, especially in the eye (danger of retinal detachment)
- Metals especially those containing ferromagnetic foreign matter
- Transdermal or other similar implants (for example, body piercings, tattoos containing ferromagnetic material as well as magnetic piercings)

Exceptions for implants, heart valves and aneurysm clips as listed above:

- If these are classified as "MR Safe", then they are not contraindicated.
- If these are classified as "MR Conditional", then special conditions apply, which are listed under exceptions below. The adherence to them is in the responsibility of the MR operator and the qualified physician.

In general, MR examinations performed despite these contraindications are the responsibility of the physician.

WARNING

Electronic and/or electrically conductive implants and magnetizable inclusions in static and low-frequency magnetic fields and RF fields!

Risk of patient death/patient injury

- ◆ Ask the patient about implants and inclusions.

- ◆ Do not perform MR examinations on patients with electronic or electrically conductive implants and magnetizable inclusions.
- ◆ Ensure that patients wearing such implants and/or inclusions remain outside the exclusion zone (0.5 mT line).

Exceptions: Certain implantable medical devices have been cleared, approved and/or licensed by the Competent Governmental Authorities and/or labeled by the device manufacturer as “MR Conditional”. For such implantable medical devices, the previously mentioned list of general contraindications and the warning may not be applicable in its entirety.

It is the responsibility of the device manufacturer to declare an implantable medical device as MR Conditional if appropriate and to define the conditions (constraints) for safe MR scanning. The MR operator must be aware of any such conditions for MR scanning. It is the obligation of the MR operator to assure that these conditions are strictly adhered to. To obtain these specific conditions the MR operator may refer to the labeling of the implantable medical device or contact the device manufacturer. Siemens Healthineers does not assume responsibility or liability for the operation of the MR system with any implantable medical device. Especially, Siemens Healthineers is not responsible for controlling technical parameters of the MR system other than those defined by the normal operating mode, the first level controlled operating mode and the data provided in the system owner manual, such as spatial gradient field plots.

CAUTION

Eddy currents induced by low-frequency magnetic fields!

Patient burns

- ◆ Do not examine patients with electrically conducting implants or prostheses.

CAUTION

Electrically conducting objects!

Injury to patient due to warming

Incorrect diagnosis due to artifacts


- ◆ Request that the patient removes all electrically conducting objects, e.g. necklaces, rings, braces, rubber bands for long hair, piercings as well as jewelry.
- ◆ Request that the patient removes all clothing including electrically conducting material, for example, bras, metallic appliqués or woven metallic yarns.
- ◆ Inform patients that eyeliners and tattoos may contain ingredients causing artifacts or skin irritations during MR examinations. In some cases, patients have been burned.
- ◆ To prevent injuries, instruct patients to remove makeup prior to the examination.
- ◆ Instruct patients to seek medical attention in case of discomfort during or following the MR examination.

Possible side effects

Possible undesirable side effects for MR are dizziness, heating, claustrophobia and nerve stimulation.

Physiological effects

Due to the presence of alternating electromagnetic fields, patients may experience various physiological effects during MR measurements:

- SAR: Warming of body tissue due to RF fields of the RF transmitter coil

These physiological effects can be quantified as the technical value SAR.

It is generally accepted that no published evidence supporting the occurrence of cumulative and/or long-term effects after exposure to EMF emitted by the MR equipment exists.

Ambient conditions

	Control room	Equipment room	Examination room
Temperature	15–30°C (59–86°F)	15–30°C (59–86°F)	18–22°C (64–72°F)
Relative air humidity	40 to 60%	40 to 80%	40 to 60%
Absolute air humidity	<11.0 g/kg	<11.0 g/kg	<11.0 g/kg
Air pressure	1060–700 hPa*		

* 700 hPa is equivalent to a height of about 3000 m above sea level

Authorized operating personnel

The MR system must be operated according to the intended use and only by qualified persons with the necessary knowledge in accordance with countryspecific regulations, for example, physicians, trained radiological technicians or technologists, having received the necessary user training.

This user training must include basics in MR technology as well as safe handling of MR systems. The user must be familiar with potential hazard and safety guidelines and emergency and rescue scenarios. In addition, the user has to have read and understood the contents of the operator manual.

Please contact Siemens Service for more information on available training options and suggested duration and frequency of such training.

Definition of different persons

Term used	Explanation
User/Operator/Operating personnel	Person who operates the system or software, takes care of the patient or reads images Typically physicians, trained radiological technicians, or technologists

System owner	Person who is responsible for the MR environment. This includes legal requirements, emergency plans, employee information and qualifications, as well as maintenance/repair.
MR Worker	Person who works within the controlled access area or MR environment User/Operator as well as further personnel (for example, cleaning staff, facility manager, service personnel)
Siemens Healthineers Customer Service/Siemens Healthineers Service/Siemens Service/service personnel	Group of specially trained persons who are authorized by Siemens Healthineers to perform certain maintenance activities References to "Siemens Service", "Siemens Healthineers Service", or "Siemens Healthineers Customer Service" include service personnel authorized by Siemens Healthineers.

4 Classification of the Medical Device

Type of the medical device in accordance with the nomenclature for classification

4.1.1 Europe (CE)

Kind of product:	Active medical device for diagnosis
Classification rules according to Annex VIII Medical Device Regulation (EU) 2017/745:	Rule 10 and Rule 11
Classification according to Annex VIII Medical Device Regulation (EU) 2017/745:	Ila
Classification criteria:	Rule 10: Active devices intended for diagnosis are classified as Class Ila if they are intended to supply energy which will be absorbed by the human body, except for devices used to illuminate the patient's body, in the visible spectrum. Rule 11: Software intended to provide information which is used to take decisions with diagnosis or therapeutic purposes is classified as class Ila

4.1.2 Russia

Classification in accordance with Russia Nomenclature classification:	135160
--	--------

Software classification according to IEC 62304:2006

The Medical Device's Software Safety Class according to IEC 62304:2006 is **Class B**.

Brand name and version of SW: *syngo MR XA50*

Class of the Medical Device in accordance with IEC 60601-1

According to IEC 60601-1:2005, the Medical Device is the **Class I**.

IP protection class

IP protection class according to IEC 60529: **IP X0**

Applied parts

According to IEC 60601-1 Cl. 7.2.10 and 8.3 the following parts are applied parts of the following types:

Patient table	Type B applied part
Local coils, e.g. Head/Neck Coil Spine Coil Contour M Coil Breast BI Coil BioMatrix Contour L Coil Contour S Coil Coil Link	Type BF applied part
Body coil (integrated):	Type B applied part

Explosion protection

The MR system is not intended for operation in areas prone to explosion (for example, highly flammable mixtures of anaesthesia gases with air or oxygen or nitrous oxide)

Cleaning system components

 **WARNING**

Flammable cleaning or disinfection agents may cause fire or explosion!

Injury to patient

- ◆ Disinfect the system components with commercial disinfecting materials. However, do not use flammable solvents such as undiluted alcohol, benzine, or acetone. Follow the manufacturer's instructions.

The following table lists classes of active agents that have been tested and are approved. The detergents listed under "not approved" must not be used.

Approved	Not approved
• Aldehydes • Quaternary compounds • Guanidine derivatives • Peroxide compounds • Pyridine derivatives • Alcohol solution (for example, 70% isopropyl alcohol and 30% water) • Chloro derivatives • Commercially available cleaning agents, detergent substances • Alkylamine • Organic acids	• Benzine undiluted • Acetone • Phenols

Clean the system and its components with a damp cloth, if not specified otherwise in the individual subchapters.

Degree of safety in the presence of a flammable anesthetics mixture with air or with oxygen or with nitrous oxide in accordance with IEC 60601-1

No AP or APG category equipment.

Operating mode in accordance with IEC 60601-1

Continuous operation with short-term load.

Mains operated equipment with additional power sources:

None

5 Technical Data: Specifications and Characteristics

Design data

The MAGNETOM Free.Max System is available as the following:

- MAGNETOM Free.Max with software *syngo* MR XA50

Photos of the device

Front view:
with standard patient table

with patient table with vertical drive

*with high load patient table**



** The outside appearance of high load patient table is the same as a patient table with vertical drive*



The list of components and accessories

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

				
8.	Пульт аварийного выключения системы.	Magnet Stop switch	A quench can be released with the Magnet Stop switch. This switch is installed on-site as well.	
9.	Выключатели аварийные(не более 5 шт., при необходимости).	Power-Off (Emergency Shutdown) switch	The voltage to the MR system can be turned off via a Power-Off (Emergency Shutdown) switch installed on-site. The switch can be used to switch off the electric power of the	

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





				
11.	Клавиатура специализированная.	PC Keyboard US English	The MR system is equipped with an original Siemens keyboard. This keyboard is a modified Windows keyboard where the numeric keys have been replaced with symbol keys.	
12.	Катушка для головы и/или шеи Head/Neck.	Head/Neck Coil	With the Head/Neck Coil the following application areas can be covered: - Head - Neck - MR Head Angiography - MR Neck Angiography - Combined head / neck	

13.	Катушка для отделов позвоночника Spine.	Spine Coil	Respiratory sensors are integrated in the BioMatrix Spine coil and measure the patient's respiratory signal in head-first and feet-first position. The sensor loops measure the change in impedance resulting from the shift of the tissue and organs during the inhaled and exhaled phase of the patient's respiration as soon as the patient is lying on the table. The BioMatrix Respiratory Sensors can be used to trigger MR sequences based on the respiratory cycles of the patient without the need and workflow impediments of a respiratory belt.	
14.	Катушка гибкая многофункциональная Contour M или BioMatrix Contour L	Contour M Coil Contour M -> L (incl. Coil Link)	The Flex coil is available in sizes small and large. Depending on the area to be examined, the coil is positioned on the patient and is used to image the area of the body it covers.	

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
	сканирования (не более 100 шт., при необходимости).		DIXON & Fast DIXON for Fatsat WARP and Advanced WARP (SEMAC) integrates different techniques tailored to reduce susceptibility artifacts caused by orthopedic MR-conditional metal implants. License key for activating software modes for collecting and evaluating volumetric data: 3D TOF 3D SPACE 3D CISS 3D contrast enhanced MRA Ultra-fast high resolution 3D pulse sequences 3D MRCP 3D VIBE 3D MEDIC 3D FLASH 3D HASTE 3D Turbo Spin Echo 3D TurboFLASH - 3D MPRAGE 3D GRE 3D FISP 3D Phase Contrast Angiography 3D TurboFLASH 3D post-processing MPR, MIP, MinIP, VRT	
--	---	--	--	--

			License key for activating software modes for viewing and analyzing dynamic series of Images: Dynamic ce-MRA for 3D imaging over time. Dynamic volume examinations with 3D VIBE Dynamic (VIBE and Quick-FatSat) License key for activating software modes for working in DICOM format: DICOM worklist support (incl. (DICOM – Export, Import, Transfer, Record to media) • DICOM send/receive • DICOM query/retrieve • DICOM SC storage commitment • DICOM basic print • DICOM modality worklist • DICOM MPPS modality performed procedure steps • DICOM structured reports • DICOM study split License key for activating program modes for transferring ready-made scanning programs: PhoenixZIP - Exchange of pulse sequences data (e.g. via Internet) by drag & drop of	
--	--	--	--	--



			clinical images. PhoenixZIP allows transfer of whole measurement programs.	
19.	Лицензионные ключи для активации дополнительных программных режимов для планирования и управления исследованием (не более 100 шт., при необходимости).	License keys to activate additional software modes for study planning and management	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. MR protocols module - Connect The MR protocols module offers one central workspace for protocol configuration and management for a fleet of scanners. MR protocols module installation The MR protocols module offers one central workspace for protocol configuration and	


			management for a fleet of scanners. Protocol modifications must no longer be made manually at each individual scanner. With the MR protocols module, protocol changes can be made independently without blocking the scanner and automatically distributed from a central workspace. Furthermore, it supports the setup of consistent examination protocols to maintain a common organ, exam, and program structure across an enterprise fleet of scanners for different clinical demands. Expert – i - This software application enables remote access to the system (connected via local area network) for planning and processing.	
--	--	--	--	--

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



				
9.	Катушки гибкие BioMatrix Contour L (не более 4 шт.).	BioMatrix Contour L Coil	The BioMatrix Contour L Coil is a sleek lightweight blanket-like coil for large field of view coverage. Coil electronics are miniaturized to enable ultimate flexibility and patient comfort. This coil is ideal for Body and Thorax imaging as well as bi-lateral hip examinations The BioMatrix Contour L hardware is ready equipped with a Respiratory Sensor. The	


			respiratory triggering functionality is not available with syngo MR XA50A. They measure the patient's respiratory signal in head-first and feet-first position. The sensor loops measure the change in impedance resulting from the shift of the tissue and organs during the inhaled and exhaled phase of the patient's respiration as soon as the patient is lying on the table. The BioMatrix Respiratory Sensors can be used to trigger MR sequences based on the respiratory cycles of the patient without the need and workflow impediments of a respiratory belt.	
10.	Катушки гибкие Contour M (не более 4 шт.).	Contour M Coil	The Contour M coil is lightweight and super flexible to be able to comfortably adapt to the patient anatomy. With minimal surfaces and interfaces, the Contour M is hygiene-friendly and easy to clean	





11.	Катушки гибкие Contour S (не более 4 шт.).	Contour S Coil	The Contour S Coil is lightweight and super flexible to be able to comfortably adapt to the patient anatomy. It provides high image quality for body regions such as knee, shoulder, hand/wrist, ankle and forefoot. With minimal surfaces and interfaces, the Contour S is hygiene friendly and easy to clean	
12.	Переходники для подключения гибких катушек (не более 5 шт.).	Coil Link	Interchangeable coil cable connecting Contour S and Biomatrix Contour L Coils to MAGNETOM Free. Platform.	
13.	Камеры для видеонаблюдения за пациентом (не более 10 шт.).	Patient Observation Camera Patient Observation Kit Desk Patient Observation Kit Wall-mount Patient Observation Camera Only	This feature provides a connection from the operator to the patient. It improves the patient experience by reducing anxiety through virtual handholding. It also allows the operator to see the patient in the scanner bore and to observe if breath commands are being followed.	

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









20.	Наушники многофункциональные, немагнитные (не более 3 шт.)	In-Ear Headset	In-Ear Headphone for easy communication with patient while using the head coil.	
21.	Роутер для получения изображений в режиме реального времени	Access-i	Access-i is an interface and Dot AddIn enabling a connection of third-party devices to the syngo Acquisition Workplace in order to receive images in real time, as well as to control the MR exam measurements. This package includes the Access-i license and a network router required for the connection of third-party workstations.	

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

25.	Лицензионные ключи для активации программных режимов для автоматизированных исследований (не более 100 шт.).	Optional license keys for activating software modes for automated research	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. myExam Brain Assist is intended to simplify general brain examinations. It allows the user to attach guidance information to measurement steps (Guidance view), to arrange most important pulse sequence parameters on one dialog (Parameters view), and to enter patient-specific information in the Patient View dialog. By selecting a predefined exam strategy, the appropriate pulse sequences will be loaded into the measurement queue. myExam Spine Assist - simplifies examinations of the spine (morphology only). It	
-----	--	--	--	--


			allows users to attach own guidance information to measurement steps (Guidance view), to arrange their most important pulse sequence parameters on one dialog (Parameters view) and to label spinal structures. myExam Large Joint Assist and was extended to also support hip and shoulder examinations and collectively are called the myExam Large Joint Assist. The myExam Large Joint Assist includes workflows for the knee, hip, and shoulder for 2D imaging techniques. Additionally, workflows for single and bilateral hip examinations are included. The AutoAlign functionality which provides prepositioning and pre-alignment of the slice groups by the system, is also available for hip and shoulder workflows. myExam Abdomen Assist is intended to simplify examinations of the liver, biliary, and pancreatic system	
--	--	--	---	--

			and kidneys (morphology only). It allows the user to attach his guidance information to measurement steps (Guidance view), to arrange his most important measurement parameters on one dialog (Parameters view), and to enter patient-specific information in the Patient View dialog. By selecting his predefined exam strategy, the appropriate pulse sequences will be loaded into the measurement queue. myExam Brain Autopilot helps users to automate intelligently. It enables less trained staff to scan brain MRI at high quality. By using automation, it takes away burdensome routine tasks for all technologists. myExam Spine Autopilot - helps users to automate intelligently. It enables less trained staff to scan general cervical, thoracic and lumbar spine MRI at high quality. By using automation, it takes away burdensome routine tasks for all technologists.	
--	--	--	---	--

			myExam Knee Autopilot helps users to automate intelligently. It enables less trained staff to scan knee MRI at high quality. By using automation, it takes away burdensome routine tasks for all technologists. myExam Shoulder Autopilot helps users to automate intelligently. It enables less trained staff to scan shoulder MRI at high quality. By using automation, it takes away burdensome routine tasks for all technologists. myExam Hip Autopilot helps users to automate intelligently. It enables less trained staff to scan hip MRI at high quality. By using automation, it takes away burdensome routine tasks for all technologists.	
--	--	--	--	--

26.	Лицензионные ключи для активации дополнительных программных режимов для экспертного сбора и обработки данных (не более 100 шт.).	Optional Application Packages;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. Neuro Perfusion Package helps to streamline the clinical workflow by inline post-processing in dynamic susceptibility contrast (DSC) based perfusion imaging. This makes it possible to see perfusion maps immediately. SWI-EPI Sequence is a high-resolution 3D imaging technique for the brain with ultra-high sensitivity for microscopic magnetic field inhomogeneities caused by deoxygenated blood, products of blood decomposition and microscopic iron deposits. Among other things, the	
-----	--	--------------------------------	--	--

			method allows for highly sensitive detection of cerebral hemorrhages and high-resolution display of venous cerebral blood vessels. The new EPI-based approach offers faster SWI with improved SNR efficiency.	
27.	Лицензионные ключи для активации дополнительных программных режимов для экспертной обработки изображений (не более 100 шт.).	Optional Postprocessing Packages.	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. syngo.MR Neuro Perfusion Engine extends the MR Neurology workflow with a complete package for advanced processing and evaluation of brain perfusion datasets syngo.MR Neuro Perfusion Mism. extends the assessment of brain perfusion datasets with	


			calculation of mismatch area and mismatch ratio between two different contrasts (for example diffusion and perfusion) to support therapy decisions.	
28.	Лицензионные ключи для активации дополнительных программных режимов сбора данных (не более 100 шт.).	Additional optional single Applications;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. StarVIBE enables the examination during continuous shallow breathing. StarVIBE can also be used in other body regions where motion is difficult to avoid.	 3D ASL acquisition output: (left) Original EPI image series; (right) Perfusion-weighted image.

29.	Лицензионные ключи для активации дополнительных программных режимов для просмотра и анализа динамических серий изображений (не более 100 шт.).	Additional optional single Applications;	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. TWIST is a sequence for time-resolved (4D) angiographic imaging with high spatial and temporal resolution, in particular for angiography. For fat suppression, water-selective excitation is used.	
-----	--	--	--	---

30.	Лицензионные ключи для активации дополнительных программных режимов для планирования и управления исследованием (не более 100 шт.).	MR protocols module - Connect MR protocols module installation	The optional functionality covers software features and applications which can additionally be bought by customers. They comprise the standard functionality and the optional licenses, depended on the individual customer order. For new MR system delivery, the software licenses are pre-installed in the MR system host computer. MR protocols module - Connect offers one central workspace for protocol configuration and management for a fleet of scanners. MR protocols module installation offers one central workspace for protocol configuration and management for a fleet of scanners.	
31.	Источники бесперебойного электроснабжения (не более 10 шт.).	UPS system UPS Battery module	UPS system Liebert GXT5 3000IRT2EXLE for MAGNETOM Free.Max and Free.Star for safeguarding computers. Including Power Cable of 9 m for connecting the UPS.	

			- Power output: 3.0 kVA / 2.7 kW - Bridge time: 5 min full load / 13 min half load - Input voltage: 230 VAC	
32.	Чиллер.	Siemens Outdoor Cooling System	Siemens Outdoor Cooling System: - Consists of two units for ultimate reliability 16kW cooling power @ 50Hz 17 kW cooling power @ 60 Hz Power consumption = 17.6 kW @ 50Hz Power consumption = 21 kW @ 60Hz Max Current Consumption (single unit) for circuit breaker selection:	


			- 3 PE ~ 50 Hz, 400 V: 15.2 A 3 PE ~ 60 Hz, 460 V: 16.1 A	
33.	Резервный теплообменник	Backup Heat Exchanger - Passive	This optional heat exchanger is for use with a Passive Cooling Set up to enable connection to a backup Chiller for the Magnet Refrigeration. This option is important if there is no reliable chilled water supply or existing redundant chiller set up at the customer site to ensure cooling to the Magnet Refrigeration.	


34.	Электропитание распределительное.	Distribution board	Distribution board	
35.	Кабина радиочастотная	High-V RF Cabin	Steel-full-room RF cabin, in a module method of building Maximum outside dimensions: - Length x Width up to 35 sqm - Height up to 3,00 m Screen damping value: 90 dB (Co-Siting 100 dB) at 15 until 65 MHz	
36.	Фильтр радиочастотной кабины.	Filter plate for RF-Cabin	Filter plate in low leakage circuitry for shielded rooms (MRI) 6 separate lines, 250 V, 50/60 Hz, 250Vdc, 32A	

37.	Облицовка радиочастотной кабины.	High-V RF Cabin Interior	Standard interior décor with LED lighting for the RF cabin. Maximum outside dimensions: - Length x Width up to 35 sqm - Height up to 3,00 m	
38.	Звукоизоляция радиочастотной кабины.	RF Cabin Sound Proofing	Basic sound protection system. Steel modular cabin (Art. No. 14418680) Standard outside dimensions of cabin: L x W x H = 7,090 x 4,090 x 2,750 mm.	
39.	Шкаф с сетевыми фильтрами.	Filter Cabinet, big	Cabinet for filter plate inside of the RF cabin. Wall-mounted cabinet custom-made. Outside dimensions appr. 1.300 x 420 x 2.700 mm (L x W x H)	

40.	Люки аварийные для экстренной эвакуации (не более 5 шт.).	Emergency hatch	With decorative plate on both sides, white, to be opened from the outside, approx. 600 x 600 mm.	
41.	Окна для радиочастотной кабины (не более 5 шт.).	Additional Day-Light Window	Observation / daylight window that can be opened, installation sideways. Dimensions: 1,200 x 800 mm.	

42.	Двери для радиочастотной кабины (не более 5 шт.).	Threshold-free door	Soundproof threshold-free door for the RF-Cabins of MAGNETOM 1.5T/3T that can be opened and closed without effort by pneumatic mechanism. Clear opening: W x H 1,200 x 2,050 mm	
43.	Стол оператора (не более 30 шт.).	Control Room Desk 1.35m	Workplace for operating personnel made of wooden, plastics and aluminum parts. Includes mount for electronic components. Length: 1350 mm Width: 800 mm Height: 720 mm Weight: 40 kg Load capacity: 100 kg	

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Technical Requirements

Technical Requirements based on the key parameters and values as presented in the chapter below.

Key parameters & values provided by manufacturer

5.1.1 Magnet system

Magnet parameters	
Operating field strength	0.55 T
Magnet type	Superconductor
Field stability over time	<0.1 ppm/h
Weight (with cryogen) ^{a)}	1550 kg
Magnet length	1.48 m
System length cover to cover	165 cm
Bore size ^{b)}	80 cm
Type of installation	Fixed
Decay characteristics from full field to 20 mT	approx. 1.6 s

a) Normal weight without any optional equipment

b) includes gradient coil, RF body coil

Guaranteed homogeneity (based on highly accurate 24 plane plot)	
10 cm DSV	0.02 ppm
20 cm DSV	0.075 ppm
30 cm DSV	0.3 ppm
40 cm DSV	1.4 ppm
50×50×45 cm ³ DEV	5 ppm

In compliance with the German «Qualifikationsvereinbarung». Standard deviation V_{rms} (Volume root-mean square) measured with highly accurate 24 plane plot method (20 points per plane). Standard active shim with 3 linear channels (1st order). DSV = diameter spherical volume (x, y and z direction). DEV = Diameter elliptical volume.

Shimming	
Both: passive and active shimming. Passive shimming during installation.	
Standard active shim with 3 linear channels (1st order).	
3D Shim	Patient-specific automated shim
	Time to shim = approx. 6 s

Shielding		
Active Shielding (AS)	7th generation active shielding (AS) technology with counter coils	
Fringe field (axial×radial)	0.5 mT ^{a)}	4.0 m×2.5 m
	0.1 mT	6.0 m×3.8 m
External Interference Shield (E.I.S.)	Patented shielding system integrated into the magnet	
	Continuous compensation and automatic suppression of external magnetic field interferences during measurement (caused by moving ferromagnetic objects or nearby power lines)	

a) pacemaker safety limit

Magnet cooling system	
Refill interval (typical) ^{a)}	Not applicable
Boil-off rate (typical) ^{b)}	0.0 l/year
Helium Inventory ^{c)}	0.7 l liquid
Minimum helium level	Not applicable
Cryostat	Stainless steel

a) For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

b) For typical clinical use, depending on sequences and operating time with running helium compressor. The system needs to be serviced at regular interval. Undisturbed magnet cooling for 24 hours and 7 days a week.

c) Liquid helium volume while operational

Static magnetic field ^{a)}	
Maximum static magnetic field	0.66 T
Maximum spatial gradient of static magnetic field	2.3 T/m

Static magnetic field ^{a)}	
Maximum product of the static magnetic field and spatial gradient of static magnetic field	1.6 T ² /m

a) Within the 80cm diameter patient accessible bore, not bare magnet figures

5.1.2 Patient handling

The patient table is available in two different configurations:

Patient table ^{a)}		
Max. weight capacity for vertical and horizontal table movement ^{b)}		220 kg (485 lbs)
Max. patient weight capacity for vertical and horizontal table movement		200 kg (440 lbs)
Max. scan range		1400 mm
Vertical table movement (only applies to patient table with vertical drive) ^{c)}	Range	480 mm – 898 mm; ± 10 mm ^{d)}
	Max. speed	40 mm/s
Horizontal table movement	Max. range	2150 mm
	Max. speed	200 mm/s
	Repositioning accuracy ^{e)}	± 1 mm
Required distance between the front of the fully extended patient table and the front wall		500 mm

High-load patient table ^{f)}		
Max. weight capacity for vertical and horizontal table movement ^{g)}		320 kg (705 lbs)
Max. patient weight capacity for vertical and horizontal table movement		300 kg (660 lbs)
Max. scan range		1400 mm
Vertical table movement ^{h)}	Range	480 mm - 898 mm; ± 10 mm ⁱ⁾
	Max. Speed	40 mm/s
Horizontal table movement	Max. range	2150 mm
	Max. speed	200 mm/s
	Repositioning accuracy ^{j)}	± 1 mm

Required distance between the front of the fully extended patient table and the front wall	500 mm
Max. initial pulling force under emergency	240 N
Max. average pulling force under emergency	240 N

- a) Including standard patient table and patient table with vertical drive.
b) Optional: vertical table drive to lower patient to 48 cm for easy access
c) Optional: vertical table drive to lower patient to 48 cm for easy access
d) Depending on the floor conditions
e) Accuracy for repositioning from one direction
f) Patient table with vertical drive only
g) Optional: vertical table drive to lower patient to 48 cm for easy access
h) Optional: vertical table drive to lower patient to 48 cm for easy access
i) Depending on the floor conditions
j) Accuracy for repositioning from one direction

5.1.3 Hearing protection data

The A-evaluated, effective sound pressure level was measured according to NEMA MS 4-2010 (National Electrical Manufacturers Association) using the maximum gradient acoustic noise (MGAN) method.

Patient noise	
Patients require hearing protection with an SNR = 19 dB or more ^{a)}	
Noise for personnel in the examination room	
Noise measured	87.7 dB(A)
Hearing protection	SNR = 6 dB or more

a) SNR = Single Number Rating

5.1.4 Gradient system

General features
Actively shielded (AS) whole-body gradient coil system
Extremely low eddy currents
Water-cooled coil and amplifier for maximum performance
All axes force compensated

Gradient amplifier ^{a)}	
• Water-cooled, highly compact, modular design • Ultra-fast solid-state technology with very low switching losses	
Max. output voltage ^{b)}	1200 V
Max. output current ^{c)}	330 A
Max. power	0.4 MW

Gradient amplifier^{a)}

1) Values for each of the 3 gradient axes

2) Values for each of the 3 gradient axes

3) Values for each of the 3 gradient axes

Gradient performance for each axis	
Max. amplitude	26 mT/m
Min. rise time	578 μ s
Max. slew rate	45 T/m/s
Vector gradient performance (vector addition of all 3 gradient axes)	
Max. eff. amplitude	45 mT/m
Max. eff. slew rate	78 T/m/s
Gradient duty cycle	100%

5.1.5 RF technology

RF transmit technology		
Frequency stability (5 min)	$\pm 2 \times 10^{-10}$	
Frequency control	32 bits (0.015 Hz)	
Phase control	16 bits (0.006 Degrees)	
Body coil	- Integrated whole body no tune transmit/receive coil with 16 rungs, 2-port feeding - Optimized RF efficiency and signal-to-noise ratio (SNR) - Real-time feedback loop for unmatched RF stability	
Transmitter path	Feedback loop for excellent RF stabilization	
	Transmit amplitude	16 bit control 25 ns resolution
	Gain stability (after first minute)	<0.05 dB (1 s) <1.0 dB (5 min) Sinc-shaped < 0.3 dB (5 min), including body coil
Transmitter amplifier	Extremely compact, water-cooled solid state amplifier, fully integrated in the electronics cabinet.	
	Transmit amplifier bandwidth	410 kHz
	Power capability of the amplifier	12 kW
	Available power at body coil input	10.2 kW

RF receiver technology	
Number of coil elements ^{a)}	51
Number of independent receiver channels that can be used simultaneously in one single scan in one FOV, each generating an independent partial image	24
Quadrature demodulation and filtering	Digital
Receiver bandwidth	500 Hz (for each channel)
Receiver signal resolution	32 bit
ADC sampling rate	120 MHz
Preamplifier noise figure	<0.6 dB
Dynamic range at coil connector (referred to 1 Hz resolution bandwidth)	151 dB/Hz instantaneous at coil connector

a) Channels that can be connected simultaneously

5.1.6 RF coils

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions (L × W × H)	Weight
Head/Neck Coil	Receive	12 channel phased array	¹ H	• Head examination • Neck examination • Combined Head/Neck examination • Temporomandibular joints (TMJ)	495 mm × 350 mm × 334 mm	6 kg
Spine Coil	Receive	9 channel phased array	¹ H	• Spine examination • Posterior part for thorax, abdomen, pelvis, hip	651 mm × 466 mm × 67 mm	5 kg
BioMatrix Contour L Coil	Receive	6 Channel phased array	¹ H	• Thorax • Abdomen • Pelvis • Bilateral hip	416 mm × 632 mm × 65 mm	0.7 kg

Coil name	Coil type	Coil design	Intended use		General data	
			Resonant nuclei	Applications	Dimensions (L x W x H)	Weight
Contour M Coil	Receive	12 Channel phased array	¹ H	• Uni lateral Hip • Shoulder • Knee • Hand (whole) • Foot (whole)	301 mm x 598 mm x 48 mm	1.25 kg
Contour S coil	Receive	6 Channel phased array	¹ H	• Shoulder • Elbow • Wrist • Hand • Unilateral hip • Knee • Foot • Ankle	269 mm x 446 mm x 40 mm	0.37 kg
Coil Link	Receive	6 Channel interface	¹ H	Use together with the following coils: • BioMatrix • Contour L coil • Contour S coil	Length: 1512 mm	0.55 kg
Breast BI	Receive	7 Channel phased array	¹ H	• Breast	561 mm x 565 mm x 244 mm	8.17 kg

5.1.7 Sequences

Basic set of sequences provided on the scanner.

Spin echo family of sequences

- SE (Spin Echo): single
- TSE (Turbo Spin Echo)
- HASTE (Half-Fourier Acquisition with Single Shot Turbo Spin Echo)
- SPACE (Sampling Perfection with Application-optimized Contrast using different flip angle Evolutions): for fast 3D spin echo imaging

Gradient echo family of sequences

- FLASH (Fast Low Angle SHot): spoiled GRE (gradient echo)
- VIBE (Volume Interpolated Breathhold Examination): for fast 3D gradient echo imaging
- MEDIC (Multi Echo Data Image Combination)
- TurboFLASH
- FISP (Fast Imaging with Steady State Precession)
- EPI (Echo Planar Imaging): with signal type SE or FID (Free Induction Decay)
- Segmented EPI: segmented in phase or readout direction
- BEAT: for ToF (Time-of-Flight) Angiography, T2-weighted imaging of the thorax and the abdomen.

5.1.8 Computer system

syngo Acquisition Workplace		
Host computer	Number of cores	6
	Processor	Intel CFL E-2226GE
	Clock rate	3.4 GHz
	Main memory (RAM)	64 GB
	Hard disk (DICOM Standard, ISO 9660)	SSD: 480 GB ^{a)}
	CD/DVD driver	Not built in, but optionally connectable by USB
	Media drive	SDHC card reader
24" standard monitor ^{b)}	High resolution widescreen monitor Automatic backlight control for long-term brightness stability	
	Screen size (diagonal)	24"
	Horizontal frequency	74 kHz
	Vertical frequency	60 Hz
	Screen matrix	1920 pixels × 1200 pixels
24" touch monitor ^{c)}	High resolution widescreen monitor with touch functionalities Automatic backlight control for long-term brightness stability	
	Screen size (diagonal)	24"
	Horizontal frequency	74 kHz
	Vertical frequency	60 Hz
	Screen matrix	1920×1200 pixels
Wall-mount touch monitor ^{d)}	High resolution widescreen monitor with touch functionalities Automatic backlight control for long-term brightness stability	

syngo Acquisition Workplace		
	Screen size (diagonal)	24"
	Horizontal frequency	74 kHz
	Vertical frequency	60 Hz
	Screen matrix	1920 × 1200 pixels
Measurement and reconstruction system	Processor	Intel CFL E-2124G
	Clock rate	3.4 GHz
	Main memory (RAM)	32 GB
	Hard disk	SSD: 480 GB
	Hard disk for raw data	SSD: ≥ 240 GB
	Reconstruction speed	≥18064 recons per second (256 ² PeFT, full FoV) ≥ 75675 recons per second (256 × 64 PeFT, 25% FoV)

a) Using Enhanced DICOM > 2200000 images with a matrix size of 256 × 256 can be stored, when acquiring image stacks with 25 slices per stack

b) A standard monitor without calibration is not suitable for diagnostic purposes. Please consider the initial acceptance testing for image display devices and the follow-up service for constancy testing on a regular base, as offered by Siemens service.

c) Optional

d) Optional

5.1.9 Network/data coupling

The MR system provides the level of safety according IEC 60950-1 outside patient environment. All equipment connected to the system's network/data couplings must also provide minimum level of safety according IEC 60950-1.

5.1.10 Basic syngo MR viewing, post-processing, and filming functionalities

Viewing

- Image display
 - Various display layouts selectable
 - Up to 3 patients can be simultaneously active
 - Image annotation and labeling
 - Configurable image text
 - Interactive movie
- Windowing
 - Freely selectable window width and center
 - Windowing (on all, succeeding, or single images)
 - Auto-windowing for optimized contrast

- ROI/VOI evaluation
 - Parallel evaluation of multiple regions/volumes of interest
 - Circle
 - Freehand
 - Statistical evaluation:
 - Area/Volume
 - Standard deviation
 - Mean value
 - Min/max values
- Range creation
 - New reconstructions based on acquired images
 - Parallel, radial, radial sliced, and curved ranges
- Position display
 - Displays measured slice positions on localizer image and selected series
- Image manipulations
 - Reversal of gray-scale values
 - Image rotation by user-defined angle
 - Image zoom and pan

Time Curve

- Time-intensity analysis

Calculations

- Arithmetic operations on images and series:
 - Addition, subtraction, multiplication, division within and across series
 - Calculated b-value and ADC maps
 - Image filter for smoothing and sharpening images
 - Motion correction (elastic) between volumes or phases
 - 2D and 3D distortion correction

Composing/Combine

- Composing of sagittal or coronal images (for example, whole-spine exams or run-offs)
- Combine different image stacks into a single one

3D

- MPR (Multi-Planar Reconstruction): Real-time multi-planar reformatting of secondary views
- MIP (Maximum Intensity Projection): 3D reconstructions of vessels from a 3D data set, or a 2D sequential slice data set (acquired with dedicated MR Angiography sequences)
- MinIP (Minimum Intensity Projection): Similar to MIP but reconstructs the minimum intensity (for example, for Dark Blood techniques)
- VRT (Volume Rendering Technique): 3D rendering with specific adaptable opacity and shading

Printing

- Connection via DICOM Basic Print
- Selectable positioning of images onto virtual film sheet
- Selectable various film layouts
- Windowing, image zoom and pan on film sheet
- Mother-in-Child display

5.1.11 Specific absorption rates

The specific absorption rate for the RF power and the rate of change of the gradient fields are checked according to the requirements of norm IEC 60601-2-33.

5.1.12 Cooling system

Two different customer specific cooling alternatives: Passive and Active.

Passive cooling system (for connection to available cooling system)	Water consumption	>70 l/min
	Water temperature	6°C - 14°C (43°F-57°F)
	Heat dissipation to water	≤ 34 kW
Active cooling system ^{a)}	Active cooling system with two Siemens Outdoor Units (SOU) will automatic adaptation to the required cooling demands (for example, different day/night mode) to decrease energy cost.	

a) Optional.

5.1.13 Line power supply

	Values				Tolerance
	With Passive cooling				
Voltage	380 V	400 V	440 V	480 V	± 10%
Frequency	50 Hz	50 Hz	60 Hz	60 Hz	±1 Hz
Connection value	42 kVA				n.a.
	With Active cooling (optional)				
Voltage	380 V	400 V	440 V	480 V	± 10%

	Values				Tolerance
	With Passive cooling				
Frequency	50 Hz	50 Hz	60 Hz	60 Hz	±1 Hz
Connection value	64.2 kVA				n.a.

For all products, line power has to be supplied via an on-site distribution panel as detailed in the Planning Guide. Room installation has to be in compliance with VDE 0100-710.

In all countries, compliance with local and national legal regulations is required. Additionally, we strongly recommend compliance with the regulations described herein – to the extent permitted by relevant local and national laws – in order to ensure the safety of operating personnel, patients, and third parties.

Power rating plate

380 / 400 / 440 / 480 V~	
3PHASE	50/60 Hz
SHORT-TIME	kVA: ...
LONG-TIME	kVA: ...

5.1.14 Power consumption

	XJ gradient
Off	7.01 kW
Standby	7.09 kW
System ready to measure	9.07 kW
Typical examination	14.22 kW

Remarks:

- All values are typical values, applicable for 400 V/50 Hz.
- Consumption for optional separator pump not included.

The power consumption measurement is based on the COCIR methodology – MRI – Measurement of energy consumption. Many variables impact power consumption, thus there can be no guarantee that each customer will achieve the same values.

5.1.15 Dimensions

Component	Width	Depth	Height	Weight
Examination Room				
Magnet 0.55 Tesla AS (incl. Helium)	2.20 m	1.48 m	1.97 m	1550 kg
Magnet in operation, incl. gradient coil, body coil, patient table, and covers	2.29 m	4.10 m	2.06 m	3150 kg ^{a)}
Standard patient table	0.72 m	2.45 m	0.89 m	180 kg
Patient table with vertical drive	0.72 m	2.45 m	0.48 – 0.89 m	240 kg
Required min. room height clearance	/	/	2.2 m ^{b)}	/
Min. transport dimensions	2.29 m	1.65 m	1.98 m ^{c)}	/
Control Room				
Standard monitor ^{d)}	0.55 m	0.22 m	0.42-0.5 m	9.8 kg
Touch monitor	0.56 m	0.32 m	0.46 m	12.3 kg
Touch monitor for wall mounting	0.56 m	0.32 m	0.38 m	8.0 kg
Mounting bracket ^{e)}	0.56 m	0.19 m	1.92 m	18.0 kg
Equipment Room				
Electronics cabinet, incl. system control, RF system, gradient power system, host computer, image processor, cooling	1.56 m	0.65 m	1.98 m	980 kg
Heat dissipation	≤ 1.5 kW ^{f)}			

a) Based on the standard patient table

b) Finished floor to finished ceiling

c) Without transport pallet

d) Monitor height adjustable

e) With touch monitor, see above touch monitor data

f) Only ventilation might be required

5.1.16 Additional System characteristics

RF data	
Power capability of the amplifier	12 kW
Transmit bandwidth	410 kHz
Receiver bandwidth	500 kHz (for each channel)

RF transmit coils		
Body coil	Max. applied RF field B1+	26 μ T
	Max. specified B1+ rms	4.55 μ T
	Distance to isocenter, where RF transmit field is reduced by	
	3 dB	159 mm
	10 dB	264 mm

Patient comfort	
Bore size ^{a)}	80 cm
In-bore lighting	can be set at 5 different levels
In-bore ventilation	can be set at 5 different levels
In-bore intercom	including loudspeaker, microphone and earphones

a) incl. gradient coil, RF body coil

5.1.17 Software

MAGNETOM Free.Max basic software information	
Name and Version	syngo MR XA50
Latest Release	25.04.2022

Nomenclature of SW versions for new SW products of the BL MR in NUMARIS/4, NUMARIS/X and SMR (syngo.via@MR) and the (multimodality workstation) according to SHS QR1

Position	1	2	3	4	5	6	7	8	9	10
Value	V	A	N	N	A	-	A	A	N	N

Determinations:

- Pos. 1 is always occupied by "V"
- Pos. 2 designates the main version. The main version begins with A and is incremented. (major version – significant change of the software architecture or user functionality)
- Positions 3 and 4 contain the version number (*version and subversion*). Begins with 10 and will be incremented according to the amount of newly introduced features as agreed by R&D and product definition.
e.g. Enhanced update includes significant cyber security update)
- Pos. 5 designates the MR system variation and/or SW update (bugfix(es) and minor change(s) without impact to regulatory approvals).
Begins with A and is incremented, whereby the following letters are reserved: "P" for MR-PET systems "U" and "R" for UHF / MR TR and for "Research" Systems

- Positions 6 to 10 are optional and reserved for Application Packages (AP), Service Packs (SP), Service Helps (SH), Security Deliveries (SD) and Hotfixes (HF).
(Here, the 6th position is always occupied by a “—” and the 7th and 8th positions read “AP”, “SP”, “SH”, “SD” or “HF”.) (e.g. also Software patch includes ordinary cyber security patch)
- Positions 9 and 10 begin with “01” and are incremented.
- No more than 10 characters are allowed.

5.1.18 Dimensions and Weight – specific for components

Cushions and positioning aids:

Component name	Dimensions	Net weight
Table pad-head part	365*425*20 mm	0.2 kg
Table pad-body part	890*490*20 mm	0.7 kg
Soft straps	900*60*15 mm	0.03 kg
Belts	900*60*15 mm	0.03 kg
Narrow belts	30*800*3 mm 30*500*3 mm	0.05 kg
Belt 1500	2000*60*15 mm	0.1 kg
Belt straps	2000*60*15 mm	0.11 kg
Large sandbags, small sandbags.	Large:150*400 mm Small:150*300 mm	8.3 kg

Set of positioning aids:

Component name	Dimensions	Net weight
Knee wedge	555*160*395 mm	0.74 kg
Back of the head support	250*95*250 mm	0.15 kg
Large triangles	230*115*175 mm	0.1 kg
Small triangles	175*96*48 mm	0.02 kg
Wedge	234*70*175 mm	0.6 kg
Gap cushion set large	750*350*16 mm	0.3 kg
Gap cushion set small	300*100*30 mm	0.02 kg
Positioning roll	Ø 210*400 mm	0.5 kg
Lordosis cushion	240*40*350 mm	0.1 kg

Head/Neck accessories:

Component name	Dimensions	Net weight
Mirror carrier with mirror head/neck coil	191*105*109 mm	0.2 kg
Pad 32	290*190*55 mm	0.2 kg
Cushion neck	290*213*60 mm	0.1 kg
Clamp cushions head	80*80*80 mm	0.02 kg
Positioning wedges	135*60*50 mm	0.03 kg

Contour accessories:

Component name	Dimensions	Net weight
Foot/ankle positioning pads for contour coils	235*195*300 mm	0.9 kg

Cushion for extremity exams	226*100*275 mm	0.13 kg
Cushion for extremity exams	226*129*275 mm	0.25 kg

Breast BI Coil & Accessories:

Component name	Dimensions	Net weight
Support cushion Breast BI 7	l = 365 mm w = 560 mm h = 100 mm	0.106 kg
Support cushion comfort Breast	l = 565 mm w = 560 mm h = 115 mm	0.380 kg
Hip support cushion Breast BI	l = 611 mm w = 562 mm h = 253 mm	1.740 kg
Filling cushion Breast BI	l = 550 mm w = 231 mm h = 89 mm	0.121 kg
Head rest Breast BI	l = 260 mm w = 220 mm h = 172 mm	1.620 kg
Armrest	l = 520 mm w = 610 mm h = 80 mm	0.950 kg
Fixation plates (lateral)	l = 104 mm w = 164 mm h = 15 mm	0.049 kg
Fixation plates (medial)	l = 177 mm w = 122 mm h = 35 mm	0.089 kg

Frame Carrier Breast BI 7 lateral	l = 280 mm w = 120 mm h = 30 mm	0.183 kg
Frame Carrier Breast BI 7 medial	l = 280 mm w = 120 mm h = 30 mm	0.182 kg
Cranio-Caudal Fixation Unit (complete)	l = 360 mm w = 453 mm h = 140 mm	2.9 kg
Cranio-Caudal Fixation Unit (individual)	l = 130 mm w = 165 mm h = 95 mm	0.275 kg

Breast BI Coil Compression Unit

Component name	Dimensions	Net weight
CC-320-Plus Cushion Breast	l = 115 mm w = 160 mm h = 15 mm	0.023 kg

EMC Declaration

5.1.19 Environment, RF shielding and RF interferences

5.1.19.1 Environment

For the MR system, two environmental categories according to IEC 60601-1-2:2014 are defined:

- Professional healthcare facility environment

This environment covers the control room with the MR console as well as the equipment room where all hardware equipment outside of the RF-shielded examination room is installed.

- Special environment: RF-shielded room

This environment covers the examination room where the magnet and all hardware equipment inside of the RF-shielded room are installed.

The MR system in an RF-shielded room is suitable for use in professional healthcare facility environment, such as clinics and hospitals.

5.1.19.2 RF shielding

An RF-shielded room has to meet the following specifications:

0.55 T:

- frequency range: 15~128 MHz
- minimum RF shielding effectiveness: > 90 dB
- RF attenuation of the RF cabin itself (without any feedthroughs or filters): > 100 dB

It is essential that the actual shielding effectiveness and filter attenuation of the shielded location be verified to assure that they meet the minimum specification. Verification of the shielding requirements is essential prior to handover of the system. The integrity of the RF shielding is the responsibility of the customer.

Recommended test methods for measurement of RF shielding effectiveness and RF filter attenuation are given in the "Planning Guide System".

5.1.19.3 RF operation of MR system

MR systems apply RF energy for diagnosis. Operation involving RF energy is performed in the RF-shielded room.

23.6 MHz	0.55 T
----------	--------

The MR system applies IQ-modulated (without encoding any data) repeated RF pulses. This results in a total bandwidth of the modulated signal of less than ± 500 kHz around the carrier frequency.

5.1.19.4 Compatible equipment

In principle, MR systems may be used in combination with compatible equipment only. For details on compatibility, see: Operator Manual MR System.

The use of the MR system adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, the MR system and the other equipment should be observed in order to verify that they are operating normally.

Please note, HF surgical equipment is not compatible with MR systems.

5.1.19.5 Accessories and cabling

All cables, transducers, and accessories are part of the MR system and cannot be removed or replaced by the user. Unauthorized replacement of accessories and cabling may affect the EMC behavior of the MR system and cause deviation from the EMC requirements of IEC 60601-1-2:2014, and furthermore may influence the functionality of the system.

The use of accessories, transducers, and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the equipment or system as replacement parts for internal components, may result in increased electromagnetic emission or decreased electromagnetic immunity of the equipment or system and result in improper operation.

5.1.19.6 RF interference

MR systems apply and receive RF energy for diagnosis. Thus, MR systems may interfere with other equipment located nearby.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the MR system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Failure to use the MR system in the specified type of shield location could result in degradation of the performance of this equipment, interference with other equipment, or interference with radio services.

The use of RF emitting equipment, which has not been released or certified by Siemens, inside the shielded location of the MR system is not allowed unless the equipment is switched to an RF silent mode before being placed in the shielded location. In case of doubt, please check the operator manual or the guidelines of the respective third-party equipment.

5.1.20 Electromagnetic emission

According to IEC 60601-2-33:2015, the MR system outside the special environment is in compliance with emission requirements of clause 7 of IEC 60601-1-2:2014. Inside the special environment, clause 7 of IEC 60601-1-2:2014 does not apply.

Emission tests	Compliance
Radiated RF emission (150 kHz–1 GHz)	EN55011/CISPR 11 Class A/Group 2
Conducted RF emission (150 kHz–30 MHz)	EN55011/CISPR 11 Class A/Group 2

The emissions characteristics of this equipment make it suitable for use in a professional healthcare facility environment, such as clinics and hospitals (CISPR 11 class A). If it is used in a residential environment, this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or re-orienting the equipment.

5.1.21 Electromagnetic immunity

5.1.21.1 General requirements

According to IEC 60601-2-33:2015, the MR system is intended for use in the electromagnetic environment specified as follows:

- The MR system outside the special environment complies with clause 8 of IEC 60601-1-2:2014.
- All operator and patient accessible parts in the special environment are compliant with the ESD requirements in clause 8.1 of IEC 60601-1-2:2014.

The customer or the user of the MR system should assure that it is used in such an environment. Normally, the MR system is used in a professional healthcare facility environment.

Immunity tests	Test level and compliance level
Electrostatic discharge, air discharge	EN/IEC 61000-4-2 ± 2.0 kV, ± 4.0 kV, ± 8.0 kV, ± 15.0 kV
Electrostatic discharge, contact discharge	EN/IEC 61000-4-2 ± 8.0 kV
Radiated RF electromagnetic fields	EN/IEC 61000-4-3 80 MHz–2.7 GHz, 3 V/m, 1 kHz, 80 % AM
Proximity fields from RF wireless communications equipment	EN/IEC 61000-4-3 385 MHz: 27 V/m 450 MHz: 28 V/m 710 MHz/745 MHz/780 MHz: 9 V/m 810 MHz/870 MHz/930 MHz: 28 V/m 1720 MHz/1845 MHz/1970 MHz: 28 V/m 2450 MHz: 28 V/m 5240 MHz/5500 MHz/5785 MHz: 9 V/m
Burst Power supply lines	EN/IEC 61000-4-4 ± 2.0 kV, 100 kHz repetition frequency
Burst Signal input/output ports	EN/IEC 61000-4-4 ± 1.0 kV, 100 kHz repetition frequency
Surge Power supply lines (line to ground)	EN/IEC 61000-4-5 ± 0.5 kV, ± 1.0 kV, ± 2.0 kV at 0°, 90°, 180° and 270°
Surge Power supply lines (line to line)	EN/IEC 61000-4-5 ± 0.5 kV, ± 1.0 kV at 0°, 90°, 180° and 270°

Immunity tests	Test level and compliance level
Conducted disturbances, induced by RF fields	EN/IEC 61000-4-6 1 kHz, 80 % AM 3 V, 150 kHz–80 MHz 6 V, in ISM and amateur radio bands between 150 kHz and 80 MHz
Magnetic fields	EN/IEC 61000-4-8 30 A/m, 50 Hz
Voltage interruptions on power supply input lines	EN/IEC 61000-4-11 Voltage interruptions: 0% UT 250/300 cycles

These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects, and people.

5.1.21.2 In situ test

If the exemption specified in IEC 60601-1-2, section 8.6 was applied for a permanently installed large MR system, the MR system may be tested in situ with the general requirements defined above except for the radiated RF immunity over the entire frequency range from 80 MHz to 6 GHz.


The MR system has been tested for radiated RF immunity only at selected frequencies, and use nearby of emitters at other frequencies could result in improper operation.

Labeling information

5.1.22 At the electronics cabinet

Product identification label (MDR)

Name plate label

MDR = Medical Device Regulation

(1) Serial number of the system

(2) Configurable device

A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations. Those individual components may be devices in themselves.

(3) Identification as medical device

(4) Unique device identification

Combination label (MDR)

MDR = Medical Device Regulation

This label indicates that there are approved combinations of the MR system described in this manual with other products. Please contact Siemens Healthcare GmbH if you have any questions regarding permitted combinations.

Artwork of the Medical device label in Russian language

Manufacturing site: People's Republic of China

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3/PE ~ 380/400/440/480 В;

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Labels according to EN 50419:2006

Within the EU, products identified with this symbol are subject to guidelines 2002/96/EC for old electrical or electronics system and 2012/19/EU respectively. Please contact Siemens Service in case of questions about returning and disposing the MR system and/or its components and accessories.


Revision label

Revision:														Model - No.:																		
X	X	02	03	04	05	06	07	08	09	10	11	12	13	14	15	16	17	18	19	20	21	22	23	24	25	26	27	28	29	30	31	32
33	34	35	36	37	38	39	40	41	42	43	44	45	46	47	48	49	50	51	52	53	54	55	56	57	58	59	60	61	62	63	64	65
66	67	68	69	70	71	72	73	74	75	76	77	78	79	80	81	82	83	84	85	86	87	88	89	90	91	92	93	94	95	96	97	98

The current revision is marked with an X in the respective field. The example above shows revision level 01 of the respective component.

5.1.23 On the magnet cover

	Product identification label
--	------------------------------

MDR	 MDR = Medical Device Regulation (1) Serial number of the system (2) Configurable device A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations. Those individual components may be devices in themselves. (3) Identification as medical device (4) Unique device identification
Combination label	
MDR	 MDR = Medical Device Regulation This label indicates that there are approved combinations of the MR system described in this manual with other products. Please contact Siemens Healthcare GmbH if you have any questions regarding permitted combinations.

5.1.24 At the magnet bore (front and back)

Risk of injury

Mandatory action: refer to operator manual

5.1.25 EAC Label*


Each unit of the MAGNETOM Free.Max delivered to Russia, Belarus and Kazakhstan will be endorsed with the EAC label right beside the system labels.


**This symbol is applicable to the system until 11.12.2022, after this time the symbol is not applied.*

5.1.26 On the coils and interfaces

Breast BI coil	
Contour M Coil	
Coil Link	

Head/Neck Coil	 Head/Neck Coil (01) 04056869281858 GTIN (240) 11407659 MODEL (422) 156 ORIGIN CN (21) 1002 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in China
Spine Coil	 Spine Coil (01) 04056869281833 GTIN (240) 11407378 MODEL (422) 156 ORIGIN CN (21) 1008 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in China
BioMatrix Contour L Coil	 BM Contour L Coil (240) 11407657 MODEL (422) 276 ORIGIN DE (21) 1007 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany
Contour S Coil	 Contour S Coil (240) 11407176 MODEL (422) 276 ORIGIN DE (21) 1009 SERIAL IVK Siemens Healthcare GmbH, Henkestr. 127 91052 Erlangen, Germany

5.1.27 Other labels

Standard patient table	 eMeRge W/O Vertical table (01) 04056869244020 GTIN (240) 10911997 MODEL (422) 156 ORIGIN CN (21) 10101 SERIAL Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in P.R. China
Patient table with vertical drive	 Patient table W/I Vertical P015 (01) 04056869244013 GTIN (240) 10911996 MODEL (422) 156 ORIGIN CN (21) 11134 SERIAL IVK Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in P.R. China  eMeRge W/O Vertical table (01) 04056869244020 GTIN (240) 10911997 MODEL (422) 156 ORIGIN CN (21) 10101 SERIAL Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center, Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA Made in P.R. China

Power consumption label

380 / 400 / 440 / 480 V~	
3PHASE	50/60 Hz
SHORT-TIME	kVA: 47
LONG-TIME	kVA: 13

5.1.28 Description of the symbols on the labels

GTIN	Global GTIN Item Number
Model	Model
	Serial number
	Manufacturer
	Catalog Number
	Do not throw away. Take it to a recycling point.
	Product compliance with European regulations
	Authorized Representative in the European Community
	Caution! Refer to the instructions for use.
	Refer to the instructions for use.
	Class item II
	Working part type BF
Rx only	Caution: Federal law restricts this device to sale by or on the order
	Medical Devices
	Date of manufacture

	Top. Indicates correct upright position
	Keep dry
	Fragile, please handle with care
	Temperature limitation

6 Sterility

MR system device is not intended to be sterile.

7 The list of International and National Standards and Directives applied to the Medical Device

IEC 60601-1:2005 + CORR.1(2006) + CORR.2(2007) + AM1(2012) EN 60601-1:2006 + Cor.:2010 + A1:2013	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
IEC 60601-2-33:2010/A2:2015 EN 60601-2-33:2010/A2:2015	Medical electrical equipment - Part 2-33: Particular requirements for the basic safety and essential performance of magnetic resonance equipment for medical diagnosis
IEC 60601-1-2:2014 EN 60601-1-2:2015	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral standard: Electromagnetic disturbances - Requirements and tests
IEC 60601-1-6:2010 EN 60601-1-6:2010	Medical electrical equipment – Part 1-6: General requirements for safety – Collateral standard: Usability
IEC 62304:2006/ AMD 1 2015 EN 62304:2006/ AC: 2008	Medical device software - Software lifecycle processes
IEC 62366-1: 2015 EN 62366: 2008	Medical devices - Application of usability engineering to medical devices
EN 1041:2008+A1:2013	Information supplied by the manufacturer of medical devices
ISO 15223-1:2016 EN ISO 15223-1:2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements
ISO 10993-1:2018 EN ISO 10993-1:2009	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process
ISO 14971:2019 EN ISO 14971:2019	Medical devices - Application of risk management to medical devices
ISO 13485:2016 EN ISO 13485:2016	Medical devices - Quality management systems - Requirements for regulatory purposes
EN IEC 63000:2018	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances
National standards	
GOST R 50444-2020	Medical instruments, apparatus and equipment. General technical requirements
GOST R IEC 60601-1-2010	Electrical medical products. Part 1. General safety requirements taking into account the main functional characteristics
GOST R IEC 60601-1-6-2014	Electrical medical products. Part 1-6. General safety requirements, taking into account the main functional characteristics. An additional standard. Operational suitability
GOST R ISO/IEC 9126-93	Information technology. Evaluation of software products. Quality characteristics and guidelines for their use
GOST 28195-89	Assessment of the quality of software facilities. General Provisions
GOST R ISO/IEC 12119-2000	Information technology (IT). Software packages. Quality Requirements & Testing
GOST R IEC 62304-2013	Medical products. Software. Lifecycle Processes
GOST R IEC 60601-2-33-2013	Medical electrical products. Part 2-33. Particular safety requirements, taking into account the main functional characteristics, for medical diagnostic equipment operating on the basis of magnetic resonance
GOST R IEC 62366-1-2021	Medical products. Part 1: Design of medical devices for serviceability


GOST R ISO 9127-94	Information processing systems. User documentation and packaging information for consumer software packages
GOST R IEC 60601-1-2-2014	Medical electrical products. Part 1-2. General safety requirements with respect to essential functional characteristics - Parallel standard - Electromagnetic compatibility - Requirements and tests
GOST R 51188-98	Data protection. Testing software for the presence of computer viruses. Model Guide
GOST R ISO 15223-1-2020	Medical products. Symbols used in the labeling of medical devices, on labels and in accompanying documentation. Part 1. Basic requirements
GOST R 52770-2016	Medical products. Safety requirements. Methods of sanitary-chemical and toxicological tests
GOST 31209-2003	Containers for blood and its components. Chemical and biological safety requirements and test methods
GOST 4011-72	Drinking water. Methods for measuring the mass concentration of total iron. item 2
GOST 31870-2012	Drinking water. Determination of the content of elements by atomic spectrometry, p. 4
MI 4.1.3166-14	Gas chromatographic determination of hexane, heptane, acetaldehyde, acetone, methyl acetate, ethyl acetate, methanol, isopropanol, acrylonitrile, n-propanol, n-propyl acetate, butyl acetate, isobutanol, n-butanol, benzene, toluene, ethylbenzene, m-, o- and p- xylenes, isopropylbenzene, styrene, α -methylstyrene in water and aqueous extracts from materials of various compositions
MI 4.1.3169-14	Gas chromatographic determination of dimethyl phthalate, dimethyl terephthalate, diethyl phthalate, dibutyl phthalate, butyl benzyl phthalate, bis(2-ethylhexyl) phthalate and dioctyl phthalate in water and aqueous extracts from materials of various compositions
GOST R 55227-2012	Water. Methods for determining the content of formaldehyde, p. 6
MVI MN 1924-2003	Method of gas chromatographic determination of phenol and epichlorohydrin in model media simulating food products
MR 1436-76	Guidelines for the determination of diphenylolpropane, as well as some phenols in its presence, during sanitary-chemical studies of products made of polymeric materials intended for contact with food products
MR 2413-81	Guidelines for the determination of epichlorohydrin in extracted extracts from polymeric materials
MR 1941-78	Guidelines for the determination of vinyl chloride in polyvinyl chloride and polymeric materials based on it, in model environments simulating food products, in food products
I-880-71	Instructions for the sanitary-chemical examination of products made from polymeric and other synthetic materials intended for contact with food products
MR 1328-75	Guidelines for the determination of caprolactam in water, air and biological media
GOST 31214-2016	Medical products. Requirements for samples and documentation provided for toxicological, sanitary-chemical tests, tests for sterility and pyrogenicity
GOST R 52770-2016	Medical products. Safety requirements. Methods of sanitary-chemical and toxicological tests
GOST ISO 10993-1-2021	Medical products. Evaluation of the biological effect of medical devices. Part 1. Evaluation and research in the risk management process

GOST R ISO 10993-2-2009	Medical products. Evaluation of the biological effect of medical devices. Part 2. Requirements for the treatment of animals
GOST ISO 10993-5-2011	Medical products. Evaluation of the biological effect of medical devices. Part 5 Cytotoxicity testing: in vitro methods
GOST ISO 10993-10-2011	Medical products. Evaluation of the biological effect of medical devices. Part 10: Irritant and sensitizing studies
GOST ISO 10993-12-2015	Medical products. Evaluation of the biological effect of medical devices. Part 12: Sample preparation and control samples



8 The list of equipment and measurement tools necessary for quality control of the medical device

The following accessories for quality assurance are available at the customer. The phantoms are not calibrated.

Measurement phantoms

	Plastic bottle, 15ml Solutes per 1000 g H ₂ O dist.: 3g MnCl ₂ * 4 H ₂ O + 5g NaCl
	Plastic bottle, 5300ml Solutes per 1000 g H ₂ O dist.: 3.75g NiSO ₄ * 6 H ₂ O + 5g NaCl
	Spherical phantom D240 Diameter: 240 mm Contents: 7300 ml Solutes per 1000 g H ₂ O dist.: 1.25g NiSO ₄ * 6 H ₂ O

Phantom holders

			
Phantom holder 27	Phantom holder 52	Phantom holder 53	Phantom holder 53
			
Phantom holder 54			

9 Transportation requirements and conditions

- The magnet can be shipped and transported worldwide with or without liquid helium. Maximum possible helium level during shipping: 95%.
- The magnet is classified as “Dangerous goods” as defined by the ADR, IMDG and IATA codes and their requirements. Liquid helium is classified under UN1963 (Helium, refrigerated liquid).
- Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.

Mechanical transport conditions

- Using an air-ride truck with or without a shipping pallet the following maximum peak acceleration measured at the magnet mounting points during transportation shall not exceed:
 - 5.0g vertical shock load applied to the cryostat mounting points.
 - 2.0g horizontal shock load applied to the cryostat mounting points in any direction.
 - Shock is defined as a dynamic force with duration of no more than 0.5 seconds.
 Transportation of the system containing liquid helium must be in accordance with the appropriate national and international codes.

Road Transport

- Road transport is allowed up to 15.000 miles.
- Road transport may be with or without an approved shipping pallet. With NO shipping pallet, the magnet must be supported using transport beams.
- Road vehicles must be fitted with air ride suspension, whether or not a shipping pallet is used.
- Always transport a magnet using vehicles fitted with air ride suspension.
- Road vehicles must be suitable for the intended usage, be properly maintained at all times, and comply with all local laws and regulations.
- Relevant requirements of the ADR code shall be satisfied.

- The vehicle and operator must be approved by Siemens.
- This specification does not constitute acceptance for use as part of a mobile system. Mobile magnets require a specific trailer design for continuous road transportation.

Air transport

- Air transport is allowed up to 30.000 miles.
- If air transport is used for any section of the route, an approved shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IATA code must be satisfied.


Sea shipment

- Sea transport is allowed up to 30.000 miles.
- If sea shipment is used for any section of the route, an approved (magnet system related) shipping pallet must be used for the entire transportation process.
- Relevant requirements of the IMDG code must be satisfied.


Rail Transportation


- Rail Transportation is not allowed!
- Transportation through the channel tunnel (UK to France) is permissible if the system is mounted on an approved air ride truck. The air ride suspension may be inactive.

Magnet specifications

	Atmospheric pressure	250 hPa – 1060 hPa (appropriate transportation kit)
	Ambient temperature	-20°C – 60°C (-4 °F-140 °F)
	Relative humidity	10%–90% non-condensing (must be kept dry if unpacked)

System components specifications

	Atmospheric pressure	250 hPa – 1060 hPa (appropriate transportation kit)
---	----------------------	--

	Ambient temperature	-20°C – 60°C (-4 °F-140 °F)
	Relative humidity	10%–90% non-condensing (must be kept dry if unpacked)

Lifting and Handling

- Lifting the magnet must always be conducted by personnel who are trained in handling components and medical equipment of this size and weight.
- Lifting during transportation including the final movement onto the installed site must be performed using a suitable sized crane or forklift capable of lifting the full magnet weight.
- Suitable straps or chains capable of lifting the full magnet weight must be used. The straps or chains must be secured to the magnet lifting points using the four shackles supplied and fitted to the magnet.

10 Shelf life

Siemens considers the subject devices exempt from having a shelf life due to the fact that the devices are manufactured upon order and is therefore not susceptible to degradation that would lead to functional failure.

Once an order is placed and a machine is manufactured, it is transported, installed and QC and tune up tests are completed on site to ensure proper function of the MR scanner. The MR scanners do not spend time on the “shelf” that would pose a risk to the performance of the device or its components.

Storage specifications

	Atmospheric pressure	250 hPa – 1060 hPa (appropriate transportation kit)
	Ambient temperature	-20 °C–60 °C (-4 °F–140 °F)

	Relative humidity	10 %–90 % non-condensing (must be kept dry if unpacked)
---	-------------------	--

11 Expected service life

The MAGNETOM Free.Max is designed for an expected service life of at least 10 years with operation under average conditions for which the system is expected to be safe for use (according to e.g. IEC 60601-1:2005).

Please be aware that the safety checks and their time intervals as specified in the System Owner Manual are mandatory and where appropriate, preventive measures have to be adopted or repairs performed, to ensure the safety of the system for this time.

The system is safe beyond the expected service life, if the necessary routine checks and periodic maintenance work are performed in accordance with the Maintenance Plan.

Siemens will inform the customer about an expected service end. Note that this does not affect legal obligations.

12 Procedure of recycling and disposal (Medical Device recycling or disposal data)

On disposing of the system or parts thereof, currently valid environmental legislation must be observed.

Examples of environmentally relevant components are:

- Accumulators and batteries
- Transformers
- Capacitors
- Monitor picture tubes
- Phantoms

For details contact your local customer service representative or your Siemens regional office.

System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

Classification in accordance with SanPiN 2.1.3684-21: Class A - Waste generated during / after the use of medical devices must be classified and disposed of (disposed of) as epidemiologically safe waste, close to MSW in composition, in accordance with current legislation.

13 Usage conditions

Installation requirements

To ensure correct installation, complete all preliminary steps specified.

For detailed information on preparing the room, transporting the system, and other installation requirements, consult the project manager.

The on-site project manager is responsible for all room and system-related preparations.

The following prerequisites must be met:

- The room must be fully prepared and cleaned.
- Power and water connections must be available.
- The internal resistance of the on-site power supply unit must correspond to the requirements.
- The transport route for the system must be defined and free of obstructions (regarding dimensions, access, load).
- Ensure that the magnet room is properly equipped.
- Ensure that the openings of the RF shield are sufficiently large.
- Ensure that there is sufficient space to install the magnet.
- Take account of narrow passages due to suspended ceilings, carriers, and columns.
- Determine the final location of the magnet in the magnet room.
- Ensure that the primary cold water supply is ready for operation.
- Block off the unloading zone to ensure that passersby are not harmed.
- Do not block escape routes or essential access routes to the hospital.
- The rooms must be properly equipped. The RF-shield, including the opening for inserting the magnet, must be installed. The magnet room must also be equipped with sufficient lighting. The floor must be clear so that the magnet and system components can be transported.
- The four corners of the magnet must be drawn on the floor of the magnet room to determine the magnet's location. Refer to the project plan.
- The floor from the unloading zone to the final location must be covered with metal plates (at least 5 mm thick).
- The tools required for installation must be available on site (steel rollers, hydraulic lifts, etc.).

Room requirements

Minimum total space requirement	
Examination Room	24 m ²
Equipment Room	24 m ²
Control Room	24 m ²
Minimum Examination room size (Technical Minimum; finished room size)	
Minimum technical room height / service height [cm] ^{a)}	220
Minimum room width [cm] ^{b)}	305
Recommended room height [cm] (minimum; not mandatory)	240
Minimum room length [cm] ^{c)}	536
Area [m ²]	16.4
Minimum Equipment Room size (including service area)	
Height [cm]	220
Width [cm]	206
Length [cm]	165
Area [m ²]	3.4
Minimum Operator Room size ^{d)}	

Height [cm]	210
Width [cm]	50
Length [cm]	80
Area [m ²]	0.4

a) The minimum Service height is not the minimum room height!

b) Technical minimum room width!

c) Technical minimum room length!

d) The size is based on the room without table and the wall mount touchable main display is configured.

Operating modes

To prevent health risk during MR measurements, several international organizations (for example, IEC) and various national health organizations have published guidelines and limit values (for example, IEC 60601-2-33 as the international safety standard for MRI). In compliance with country-specific approval guidelines, they are the basis for the monitoring functions integrated in the MR system with respect to stimulation and warming effects. The limits against too intense stimulation and warming effects (for example, dB/dt limits and SAR limits) are based on current scientific literature related to safety.

Two different operating modes are available depending on the patient's tolerance. With respect to stimulation and warming effects, the operating modes can be defined and selected separately.

- Normal Operating Mode

The Normal Operating Mode can be used safely for all patients. This is the standard mode. Routine patient monitoring is required.

Scanning of pregnant patients with the Body coil must be limited to the Normal Operating Mode with respect to the SAR level.

A fetus is especially sensitive to potential thermal events during the first three months of pregnancy. Therefore, avoid scanning pregnant patients in the first three months of pregnancy and avoid scanning patients with unknown pregnancy status.

- First Level Controlled Operating Mode

In the First Level Controlled Operating Mode, patients may experience noticeable stress levels depending on the measurement programs selected. The decision to change to the First Level Controlled Operating Mode must be based on a medical consideration of the potential risks and benefits for the patient.

CAUTION

Exposure to RF electromagnetic fields in the First Level Controlled Operating Mode!

General or local hyperthermia of the patient

- ◆ Do not examine patients with restricted thermoregulatory capability (e.g. small children, elderly, sick, or medicated patients).
- ◆ Do not examine patients unable to communicate potential overheating effects (e.g. small children, seriously ill, paralyzed, unconscious, sedated, or handicapped patients).
- ◆ Ensure that patients wear light clothing (e.g. light pajamas or nightgown).
- ◆ Remove all additional insulation, e.g. blankets which could interfere with heat dissipation.
- ◆ Carefully observe the patient and advise the patient once again about the squeeze bulb.

Ensure medical monitoring of the patient (as required by the IEC standard). Also consider the need for breaks during the measurements, to allow the patient to cool off, for example.


CAUTION

Stereotactic frames and similar devices: tips of screws may heat up considerably, especially if MR examinations are performed in the First Level Controlled Operating Mode!

Localized burns of the patient

- ◆ Please observe the recommendations and notes of the manufacturer of the stereotactic frame.
- ◆ If the device consists of conductive material, only perform measurements in the Normal Operating Mode.
- ◆ If you still need to switch to the First Level Controlled Operating Mode, please observe the related safety notes.

Switching operating modes: To switch from the Normal Operating Mode to the First Level Controlled Operating Mode, the user must explicitly select and confirm the change. The request appears at the syngo Acquisition Workplace. In the First Level Controlled Operating Mode, medical supervision is mandatory.

Information about maintenance (general and on the territory of Russia)

Routine checks and periodic maintenance are necessary to ensure safe and satisfactory operation of your system.

Proof of these activities is required by the authorities in certain countries.

MAGNETOM maintenance plan

This maintenance plan shows in tabular form the activities that have to be performed by qualified and authorized service engineers in the course of periodic maintenance work.

The maintenance work is subdivided into:

- Safety check (including safety-related tests)
- Preventive maintenance
- Quality and function tests

Each table contains an introductory explanation.

Detailed working instructions for all maintenance work are provided in the service documentation for this system. These documents are not included in the shipment of the system.

The RF-Cabin and the door of the RF-Cabin are not part of the MAGNETOM system maintenance and not considered in these instructions. The RF-Cabin and the door of the RF-Cabin door have to be maintained according to the instructions of the manufacturer.

Safety check (including safety-related tests)

The following checks are intended to ensure the safety of the system. Where appropriate, preventive measures have to be adopted or repairs performed. The points to be checked are generally regulated by laws and standards.

The Safety-related Tests according to §6 MPBetreibV (Germany) are mandatory.

The specified checks and intervals correspond to the minimum requirements. Compliance with stricter national legislation may occasionally be necessary.

Object or function	Reason	What is checked	Interval
General tests			
System	Safety of patient, personnel, and system	Visual inspection of the system	Annually
Visible cabling, cable routing	Safety of patient, personnel, and system	Visual inspection of cables, cable routing	Annually
Options	Safety of patient, personnel, and system	Visual inspection of options	Annually
Phantom	Safety of patient, personnel, system, and environment	Phantoms cannot have any defects or air bubbles. Replace or refill as necessary.	Annually
Accessories	Safety of patient, personnel, and system	Visual inspection of accessories ^{a)}	Annually
Magnet	Optimum function within specifications Safety of patient, personnel, and system	Check for leaks	Annually
		Ramping down the magnet	Every 2 years
Patient fan	Safety of patient and personnel	Replace filter	Annually
Operating manuals	Safety of patient, personnel, and system	Required operating manuals are available and legible	Annually

User icons, button labeling, warning labels, warning signs	Safety of patient, personnel, and system	Visual inspection of the user icons, button labeling, warning labels and warning signs • Checking the warning signs of the patient table • Checking the hearing protection sign	Annually
Door to the examination room	Safety of patient and personnel	Checking door lock	Annually
Controlled access area for magnetic field	Prevent hazards caused by magnetic field	The controlled access area (0.5 mT area) is identified The examination room is identified with the warning "strong magnetic field"	Annually
Electrical tests			
Protective earth resistance system	Protecting personnel from electrical shock. (The protective earth resistance can change during system operation, i.e. through oxidation, corrosion, and cable breaks.)	Protective earth resistance of the entire system	Every 2 years
Function test			
Emergency power off circuit	Safety of patient, personnel, and system	Checking function of the emergency power off switches	Every 2 years
Gradient supervision	Optimum function within specifications Safety of patient, personnel, and system	Function of the gradient supervision	Every 2 years
Ramp cable supervision	Optimum function within specifications Safety of patient, personnel, and system	Function of the ramp cable supervision	Every 2 years
Magnet stop function (ERDU)	Safety of patient, personnel, and system	Operation of the magnet stop buttons	Every 2 years
Patient Table	Safety of patient, personnel, and system	Checking the distance between the patient table and the wall in the RF-room Operation of the emergency stop function Operation of the emergency release Checking tolerance limits Inspect table movement and end switches Checking the patient table protective cover	Annually
Squeeze bulb	Optimum function within specifications	Function of squeeze bulb	Annually
Quality assurance measurements (General QA)	Optimum function within specifications Patient protection	Interaction of various system components according to the guaranteed system characteristics	Annually

a) According to manufacturer recommendations

Leakage Current Measurements

National regulations have to be observed, e.g. IEC 62353 / DIN EN 62353 in Germany.

13.1.1 Instructions for measuring the equipment leakage current

The leakage current measurement for the MRI system can be omitted.

Reasons:

- the protective conductor is permanently connected.
- compliance with the **on-site protective standards**, e.g. DIN VDE 0100-710, is described in the planning documents. The **system owner** has to comply with the **on-site protective standards**.
- The protective earth resistance is checked at beginning of putting the MR System into operation, at recurring maintenance and, if necessary, after repairs, upgrades and changes.

13.1.2 Instructions for measuring the patient leakage current at the patient table

Patient table (applied part of type B)

The leakage current measurement for the MRI patient table (applied part type B) can be omitted.

Reasons:

- Applied parts with non-conductive surfaces, separation between live parts and applied parts with metal parts which are directly connected to the protective conductor according to IEC 60601-1, section 8.5.2 (3rd edition).

13.1.3 Instructions for measuring the patient leakage current at the MR local coils

Local coils (applied part of type B or BF)

The leakage current measurement for local MRI coils (applied part type B or BF) can be omitted.

Reasons:

- applied part with non-conducting surface, separation between live parts and applied parts with double isolated power supplies (two MOPP) and no occurrence of hazardous voltages in the applied part according to IEC 60601-1, section 8.5.2 (3rd edition).

13.1.4 Instructions for measuring the patient leakage current at the body coil

The leakage current measurement for the body coil can be omitted.

Reasons:

- Parts with non-conductive surfaces, separation from mains connection by double isolated power supplies (two MOPP) IEC 60601-1, section 8.5.2 (3rd edition).

Service contacts

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract. It is also possible to take out a service contract based on system usage. Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubninskaja Str. 96,

115093 Moscow, Russia
Tel.: +7 495 737 12 52
www.siemens-healthineers.com/ru/

Customer's support center (24 hours):

Phone: 8 800 550 50 22
Email: service.ru@siemens-healthineers.com

Information about calibration of the MD / requirements to calibration (if applicable)

There is no calibration required for the MR systems. Nevertheless, there are different checks and tests which have to be done as described in the system owner manual.

14 Warranty

The system owner is responsible for the software licenses including software warranty conditions. The specific warranty conditions are usually defined in the contractual agreement between Siemens Healthcare GmbH and the system owner. For Russia there may be a 14 months warranty after the signature of the delivery note.

15 Contact for reclamation

In case of any questions or reclamations, the following shall be contacted:

Germany	Russia
Siemens Healthcare GmbH Henkestr. 127, 91052 Erlangen, Germany Tel.: +49 9131 84-0 www.siemens-healthineers.com	Siemens Healthcare LLC Dubininskaya Str. 96, 115093 Moscow, Russia Tel.: +7 495 737 12 52 Fax: +7 495 737 13 20 www.siemens-healthineers.com/ru/

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