More than
75
Years

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Thank you for your expression of confidence in the KARL STORZ brand. Like all of our other products, this medical device is the result of years of experience and great care in manufacture. You and your organization have decided in favor of a modern, high-quality piece of equipment from KARL STORZ.

This instruction manual enables you to properly install, connect, and operate this medical device.

- Failure to observe and follow this instruction manual and the instruction manuals of medical devices used in combination can result in injury to patients, users and third parties as well as damage to the medical device.
- Please read all relevant instruction manuals and always follow the instructions given. Check the functioning of the medical devices used in combination.
- Take note of all the instructions on the device.
- Keep this instruction manual available for ready reference in a clearly visible location near the device.

KARL STORZ reserves the right to make engineering modifications, in particular in the context of ongoing development and improvement of our devices.

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It is recommended that the suitability of the products for the intended procedure be checked prior to use.

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1	Intended use	3
1.1	Intended use.....	3
1.2	Indications	3
1.3	Contraindications.....	3
1.4	Target groups	4
1.5	User qualification	4
1.6	Patient population.....	4
1.7	Intended application areas.....	4
1.8	Intended conditions of use.....	4
1.9	User profile of physician and medical support staff	4
2	Safety	5
3	Description	7
3.1	Device overview	7
3.2	Symbols employed	8
3.2.1	Safety instructions	8
3.2.2	Pictograms/packaging symbols	9
4	Initial operation	10
4.1	Unpacking	10
4.2	Scope of supply	10
4.3	Required accessories	10
4.4	Preparing the device for use	11
5	Operation	12
5.1	Operating the user control buttons	12
5.2	Using accessories.....	13
5.2.1	Mechanical accessories.....	13
5.2.2	HF accessories.....	13
6	Reprocessing	15
6.1	Overview table	15
6.2	General safety instructions.....	16
6.3	Accessories	17
6.4	Preparation for cleaning and disinfection	17
6.5	Manual precleaning.....	17
6.6	No use of ultrasound	18
6.7	Machine cleaning and disinfection	18
6.8	Assembly, inspection and care	18

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6.9	Packaging systems.....	19
6.10	Sterilization	19
6.11	Limits of reprocessing.....	19
6.12	Care and handling of the video connecting cable	20
7	Service, repair and disposal	21
7.1	Warranty	21
7.2	Maintenance	21
7.3	Repairing device parts	21
7.4	Repair program.....	22
7.5	Infection prevention.....	22
7.6	Disposal.....	22
8	Technical description	23
8.1	Technical data.....	23
8.2	Operating/storage conditions.....	23
8.3	Directive compliance.....	23
9	Accessories.....	24
10	Electromagnetic compatibility user information	25
11	Subsidiaries.....	43

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1 Intended use

1.1 Intended use

Mediastinoscopes serve to create access for diagnostics and treatment in the middle mediastinum. Mediastinoscopes are designed for transient use in surgically invasive interventions.

1.2 Indications

- Non-small cell lung cancer/metastases
- Small cell lung cancer/metastases
- Staging of mediastinal lymph node metastases
- Interstitial lung disease
- Esophagectomy

1.3 Contraindications

There are no known contraindications relating directly to the device. The video mediastinoscope must not be used if, in the opinion of the responsible physician, the patient's health is endangered through its use, e.g., due to the patient's general condition or if the endoscopic method as such is contraindicated.

The following contraindications also apply:

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1.4 Target groups

This instruction manual is intended for physicians, medical support staff and sterilization assistants.

1.5 User qualification

Only physicians and medical support staff with a relevant specialist qualification and who have received training on the device may use the video mediastinoscope. Only technically suitable persons may provide training.

1.6 Patient population

Gender: irrelevant
Age: irrelevant
Weight: irrelevant
State of health: suitable for endoscopic interventions, in the physician's judgement

1.7 Intended application areas

Application area of mediastinoscopy, see intended use

1.8 Intended conditions of use

The video mediastinoscope may only be used in hospitals in suitable ambient conditions. For ambient conditions see the technical data.

1.9 User profile of physician and medical support staff

- Recognized medical skills in the field of application (specialist physician, qualified medical staff).
- Adequate powers of comprehension to rationally assess the surgical situation in hand.
- Adequate language skills in at least one of the languages used in the instruction manual
- Be thoroughly trained in the operation and use of the device.
- Knowledge of the contents of the instruction manual.
- No physical impairments that could diminish perception of activation and alarm signals (visual and acoustic).

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- *Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.*
- *Reprocess reusable instruments before and after each use. When doing so, use validated procedures.*

Combination of system components

Failure to observe and follow this instruction manual and the instruction manuals of products used in combination can result in injury to patients, users and third parties as well as damage to the product.

- *Please read all relevant manuals and always follow the instructions given precisely.*
- *Check the functioning of the products used in combination.*

The video mediastinoscope can only be used as a system in combination with other instruments and devices from KARL STORZ.

- *Suitable system components are listed in Section 4.3 'Required accessories'.*
- *Combination with non-suitable instruments and devices can result in uncontrolled conduct and injury to the patient, user or third parties.*

Improper use

Incorrect application of medical instruments poses a risk of injury for patients.

- *Users of medical instruments must have an appropriate medical qualification and be acquainted with the application.*

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Risks due to damaged parts

Incorrectly assembled and damaged instruments can lead to injuries to the patient.

- *Check the following properties of instruments before and after every use to ensure that they:*
 - are free from damage
 - are in full working order,
 - have no unintentional rough surfaces,
 - have no sharp corners, burred edges,
 - have no projecting parts,
 - are complete
- *Do not leave missing or broken-off components inside the patient.*

Risks from incorrect maintenance and installation

Maintenance and installation errors endanger the safety of physicians, medical support staff and patients. The instrument doesn't work reliably or works with only with limited output.

- *Only staff authorized by KARL STORZ may perform the following work:*
 - Assembly
 - Expansion
 - Readjustments
 - Modifications
 - Repairs
- *The electrical installation of the location in which the medical device is connected and operated must conform to the applicable statutes and standards*

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
- ① Video connecting cable, fixed in place
- ② User control button: for scrolling (downwards) through the menu or activating a camera function
- ③ User control button: for accessing the menu/ making a selection
- ④ User control button: for scrolling (upwards) through the menu or activating a camera function

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3.2 Symbols employed

3.2.1 Safety Instructions

HAZARD SYMBOL: This is a hazard symbol. It warns of the risk of injury. To prevent injuries, follow all the measures marked with the hazard symbol.

Signal words:

WARNING! designates a possible imminent risk. If this is not avoided, it could lead to death or serious injuries.

CAUTION! designates a possible imminent risk. If this is not avoided, it could lead to minor injuries.

ATTENTION! designates a possibly harmful situation. If this is not avoided, the medical devices could be damaged.

NOTE! indicates special information on the operation of the device or important information.

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**3.2.2 Pictograms/
packaging symbols**

	Manufacturer
	Consult instructions for use
	Follow instructions for use
	Catalogue number
	Serial Number
	Number of products in the product packaging
	Keep dry
	Applied part type BF
	This device has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).
	CE mark
	Federal (USA) law restricts this device to sale by or on the order of a physician.

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- Kompatibler HD-Monitor.

4 Initial operation

NOTE: The video mediastinoscope including the equipment connected may only be used in medical rooms in which the electrical equipment has been installed in accordance with applicable national regulations and standards.

4.1 Unpacking

Carefully remove the video mediastinoscope and the accessories from the packaging. Check for missing items and for evidence of shipping damage. If necessary: Notify the manufacturer or the supplier. TIPP: If possible, retain the original packing materials for later use; these can be used if the device has to be transported.

4.2 Scope of supply

- Video mediastinoscope 10973 HD
- Instruction manual 96086006 DER

4.3 Required accessories

ATTENTION! Non-compatible accessories can lead to malfunctioning or damage to the video mediastinoscope.

- Only use the compatible accessories designated by KARL STORZ.

The video mediastinoscope 10973 HD requires additional accessories for normal use:

- IMAGE1 S CONNECT® TC 200
- IMAGE1 S X-LINK TC 301

NOTE: These devices have their own instruction manual (Art. no. 96206286DER)

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
5 Operating instructions

! WARNING!

Risk of injury and risk of damaging the products: The highly concentrated light energy can cause the distal end, the light port, adjacent components and tissue in front of the light emission window to heat up. This can cause burns to the patient, user and operating accessories.

- > Avoid direct tissue contact with the distal end.
- > Do not rest the video mediastinoscope on the patient or in direct contact with surgical accessories.

! WARNING!

Risk of injury: Optical radiation poses a risk of injury to eyes.

- > Never look into the free end of either a connected light cable or endoscope.
- > Always wear standardized protective equipment when using lasers.

5.1 Operating the user control buttons

② Scroll down:

Scroll downwards through menus, options, etc. The button also be assigned a programmable function/feature.

③ Menu/selection:

Call up menu. Once the menu is displayed, this button can be used to select the highlighted option (identical to Enter on a computer keyboard). Hold this button down for seven seconds to call up the 'Options' menu.

④ Scroll up:

Scroll upwards through menus, options, etc. This button also be assigned a programmable function/feature.

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Further information on the programming and operation of camera head button can be found in the instruction manual for KARL STORZ IMAGE1 S™ (Art. no.: 96206286DER).

5.2 Using accessories

The spatula with an opening at the side enables the use of accessories such as forceps, scissors and HF coagulation instruments.

5.2.1 Mechanical accessories

WARNING!

Overloading by exerting too much force may cause the medical device to break, bend, and malfunction, and consequently injure the patient or user.

- Do not overload the instruments.
- Do not bend bent instruments back to their original positions.

5.2.2 HF accessories

WARNING!

Risk of injury and risk of damage to the products: If applied parts are used outside of the field of vision, there is a risk that tissue and accessories could be damaged unintentionally.

- Always hold the applied parts of the active electrode, laser and other instruments that transmit energy pointed at the target in the field of vision during use.

WARNING!

Risk of injury: Patient leakage currents may accumulate if endoscopes and endotherapy devices supplied with energy are used together. This is particularly important if an endoscope with a Type CF applied part is used. In this case, an endotherapy device with a Type CF applied part must be used in order to minimize the total patient leakage current.

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 WARNING!

Risk of injury: The electromagnetic interference fields generated during the use of HF surgery can impair the functioning of other electronic devices. In the case of patients with cardiac pacemakers or other active implants, there is the risk of disturbances or damage.

- *If in doubt, seek expert advice from KARL STORZ.*

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- This reprocessing step is a key part of the validated process
- This reprocessing step may be carried out additionally in order to aid cleaning
- X On the grounds of potential damage, this reprocessing step must NEVER be performed
No symbol: This reprocessing step has neither been validated in respect of efficacy nor in respect of material compatibility by KARL STORZ.

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6.2 General safety instructions

WARNING!

Incorrectly reprocessed instruments expose patients, users and third parties to a risk of infection as well as the danger of an instrument malfunction. The product-accompanying documents must be observed.

- > An overview of the procedures validated by KARL STORZ can be found in the manual 'Cleaning, disinfection, care and sterilization of KARL STORZ instruments' (Art. no. 96216003 DER).
- > This manual can be downloaded or requested by visiting www.karlstorz.com.

WARNING!

The instruments are not sterile when delivered. The use of non-sterile instruments poses a risk of infection for patients, users and third parties.

- > Inspect reusable instruments for visible contamination. Reprocess contaminated instruments.
- > Reprocess reusable instruments before and after each use. When doing so, use validated procedures.

WARNING!

Risk of infection: When carrying out any work on contaminated medical devices, the guidelines of the Employers' Liability Insurance Association and equivalent organizations striving to ensure personal safety as well as other country-specific laws must be observed.

WARNING!

If CJD (Creutzfeldt Jakob disease) is suspected or indicated, KARL STORZ recommends disposing of the product in accordance with national regulations and not using it again.

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ATTENTION: Using chemicals not approved by KARL STORZ may damage medical devices.

- > Only use chemicals approved by KARL STORZ for reprocessing. A list of approved chemicals can be requested from hygiene@karlstorz.com.

ATTENTION: Prolonged immersion and incorrect concentration of the cleaning/disinfection solutions may result in damage.

- > When preparing and using solutions, follow the chemical manufacturer's specifications, paying close attention to proper concentration, exposure time and service life.
- > The microbiological range of action of the chemicals used must be taken into account.

NOTE: National laws and regulations must be observed.

6.3 Accessories

Accessories required for carrying out reprocessing:

- Brush: 27652

6.4 Preparation for cleaning and disinfection

Heavy soiling, corrosive solutions and pharmaceuticals must be removed from the medical device immediately after use. To this end, preclean the medical device by wiping down and rinsing, for example. As a general rule, KARL STORZ recommends manual precleaning under cold running water.

6.5 Manual precleaning Brushing the surfaces

Any visible contamination or heavy soiling must be removed by cleaning the surfaces under cold running water with the aid of a brush or sponge.

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6.6 No use of ultrasound

For technical reasons, the medical device is not suitable for ultrasound treatment.

6.7 Machine cleaning and disinfection

NOTE: In order to protect the contacts of the video connecting cable connector, machine cleaning and disinfection is only possible if the protective cap is fitted (see Fig.).

The following methods for machine decontamination have been validated and approved subject to compliance with the process parameters described in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments':

• Machine cleaning/thermal disinfection

Thermal disinfection is preferred. The relevant national requirements and the A_0 value must be taken into account when using this method.

The selection of a suitable slide-in tray or instrument holder, which should ensure that the medical device is thoroughly rinsed out or through, must take place in consultation with the manufacturer of the device.

NOTE: If necessary, the instrument must be dried off afterwards by hand.

6.8 Assembly, inspection and care

The cleaned and disinfected medical device must be visually inspected for cleanliness, completeness, damage and dryness:

- If residues or contamination are still present, the medical device must be manually cleaned and subjected to a full cleaning and disinfection procedure once more.
- Damaged or corroded medical devices must be withdrawn from use.
- Afterwards, a functional check must be carried out.

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NOTE: Use items from the catalog 'Care, Sterilization and Storage Techniques' during care procedures.

6.9 Packaging systems

Only standardized and approved packaging materials and systems may be used (EN 868 Parts 2-10, EN ISO 11607 Parts 1 + 2, DIN 58953).

6.10 Sterilization

The procedures as well as the process-relevant parameters for the individually validated methods are described in detail in the manual 'Cleaning, Disinfection, Care and Sterilization of KARL STORZ Instruments'. The method must be selected taking into account the respective applicable national requirements and in consultation with the device manufacturers.

The following sterilization procedures have been validated and approved by KARL STORZ for this medical device:

- **Steam sterilization using the fractionated pre-vacuum procedure**

NOTE: Sterilize greased components in a disassembled state.

The medical device must be sterilized in its fully assembled state using the fractionated prevacuum procedure (DIN EN ISO 17665-1) with the following parameters:

- 132°C with a minimum exposure time of 4 minutes and maximum of 18 minutes.

NOTE: This validation is applicable for reusable tube sets which are open at both ends and with a maximum length of 200 cm.

- 134°C-137°C with a minimum exposure time of 3 minutes and maximum of 18 minutes.

6.11 Limits of reprocessing

The end of the product's service life and correct functioning are determined by mechanical stress and chemical influences within the scope of reprocessing and application.

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6.12 Care and handling of the video connecting cable

To maximize the life of the video mediascope, please follow the guidelines given below for proper handling and care.

DO NOT PLUG THE VIDEO CONNECTING CABLE INTO THE PROCESSOR RECEPTACLE IF IT IS WET. The connector must always be completely dry and free of residue. Where necessary, remove any residual moisture with a low lint disposable cloth before plugging in the unit.

NEVER PICK UP THE VIDEO MEDIASTINOSCOPE BY ITS CABLE! Always pick up the set by the instrument handle first, followed by the cable and connector. The cable contains tiny electronic conductors and should be treated gently and carefully.

PROTECT THE VIDEO MEDIASTINOSCOPE FROM DIRECT SUNLIGHT AND EXCESSIVE HEAT.

HANDLE THE SYSTEM WITH CARE, avoiding extreme impacts to the device. When transporting, use the original packaging to prevent damage to the device.

DO NOT STORE THE VIDEO-MEDIASTINOSCOPE WITH CABLE COILED TIGHTLY OR GATHERED INTO FOLDS. Store the cable in loops at least 15 cm in size.

DO NOT RAPIDLY PULL THE CABLE WHEN UNCOILING. Harmful 'kinks' can result. Always allow the cable to uncoil gently and loosely.

NOTE: We always recommend drying with sterile compressed air after manual reprocessing and, if necessary, also after machine reprocessing.

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7 Service, repair and disposal

7.1 Guarantee

The warranty conditions provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your competent KARL STORZ subsidiary (see the section entitled 'Subsidiaries'), even during the warranty period.

Opening of the equipment or performance of any repairs or modifications to the device by unauthorized persons shall relieve KARL STORZ of any liability for its performance. Any such opening, repair or modification performed during the guarantee period shall void all warranty.

7.2 Maintenance

No further repair measures are required beyond reprocessing.

The legal requirements as per the Medical Devices Operator Ordinance (MPBetreibV) for devices to be used in combination apply (see Section 4.4). The relevant warnings in the instruction manuals must be observed.

7.3 Repairing device parts

WARNING!

Danger of electric shock. Do not open the device!

- Refer servicing to the manufacturer or the manufacturer's authorized service personnel (see Section 3 of the Medical Device Operator Ordinance).
- Removal of covers by unauthorized personnel will void the equipment's warranty.

Only personnel authorized by the manufacturer may repair device parts. The operations described in this manual are exempt from this rule. The authorized personnel may only repair device parts which the manufacturer designates as repairable.

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Upon request, the manufacturer will provide authorized personnel with information on repairs, e.g., circuit diagrams, spare parts lists, setting instructions.

The fact of having such documents at one's disposal does not correspond to authorization.

7.4 Repair program

For this medical device, individual repair is necessary. To bridge the repair period, you will receive a device on loan, which must then be returned as soon as you receive the repaired medical device.

In the case of repairs, please contact your local KARL STORZ subsidiary (see section 'Subsidiaries') or authorized dealer.

In Germany, you can refer repairs direct to:

KARL STORZ SE & Co. KG
Abt. Reparaturservice
Take-off Gewerbepark 83
78579 Neuhausen

7.5 Infection prevention

To prevent infections, contaminated devices and instruments must not be returned.

Decontaminate instruments and devices before returning. KARL STORZ returns contaminated instruments to the sender.

7.6 Disposal

This unit has been marked in accordance with the European Directive on waste electrical and electronic equipment (WEEE).

At the end of its useful operating life, dispose of the device as electronic scrap.

Please ask KARL STORZ, a KARL STORZ subsidiary (see the section 'Subsidiaries'), or your authorized dealer about the collection point in your area.

Within the scope of application of this Directive, KARL STORZ is responsible for the proper disposal of this device.

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8 Technical description
8.1 Technical data

Direction of view	30°
Angular aperture	50°
Spatula working length	198 mm
Distal outer diameter of spatula	16.85 mm
Distal inner diameter of spatula	14.85 mm
Illumination	Integrated LED source
Operating temperature	5 – 30°C
Weight with cable	590 g
Dimensions max. without cable in mm (W x H x D)	198 x 144 x 25
Cable length	2.90 m

8.2 Operating/storage conditions

	Temperature	Rel. Humidity
Storage	-10°C ... +60°C	5% ... 95%
Operation	+5°C ... +30°C	0 % ... 95 %
Air pressure	700 hPa – 1060 hPa	

8.3 Directive compliance
In accordance with the Medical Device Directive (MDD):

This medical device bears a CE mark in accordance with MDD 93/42/EEC.

A code number after the CE mark indicates the responsible notified body.

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9 Accessories
MEDIAFIT forceps

MEDIAFIT Sponge and Dissecting Forceps acc. to MAASSEN, fenestrated

MEDIAFIT Biopsy Forceps, oval jaws, size 8 x 16 mm

MEDIAFIT Biopsy Forceps, oval jaws, size 6 x 12 mm

MEDIAFIT Biopsy Forceps, spoon-shaped jaws

MEDIAFIT Grasping and Dissecting Forceps, upturned jaws

MEDIAFIT Grasping and Dissecting Forceps, long and fenestrated jaws

MEDIAFIT clip forceps, jaws curved to right

MEDIAFIT clip forceps, jaws curved to left

Coagulation Suction Tube

MEDIAFIT Outer Sheath with cleaning connection, diameter 5 mm, length 21 cm

HF instrument for mediastinoscopy

HF instrument for mediastinoscopy

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10 Electromagnetic compatibility user information

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is suitable for use in the direct proximity of an active HF electrosurgical unit in professional health institutions. Professional health institutions include medical practices, dental practices, facilities for restricted care, independent surgical centers, independent birthing centers, several treatment facilities, hospitals (emergency wards, patient rooms, intensive care wards, operating theaters, outside of the HF-shielded room of an ME system for magnetic resonance imaging).

WARNING!

The use of this device next to or with other devices should be avoided as this can result in incorrect operation. If such use is necessary, this and the other devices should be monitored to ensure that they are functioning correctly.

WARNING!

The use of an accessory or cables for the video mediastinoscope 10973 HD with the IMAGE1 S™ other than those specified in this manual may result in increased emissions or decreased immunity of the video mediastinoscope 10973 HD together with the IMAGE1 S™. When using an accessory or cables with the video mediastinoscope 10973 HD together with the IMAGE1 S™ other than those specified in this manual, it becomes the responsibility of the user to determine compliance of the video mediastinoscope 10973 HD and IMAGE1 S™ with IEC 60601-1-2 when using this item.

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 WARNING!

Portable HF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the video mediastinoscope 10973 HD in conjunction with the IMAGE1 S™, including cables specified by the manufacturer. Otherwise, performance may be impaired.

NOTE: Based on the device's emissions properties, it is suitable for use in industrial areas and hospitals (CISPR 11 Class A). If the device is used in a domestic environment (for which CISPR 11 Class B is normally required), the device may not offer sufficient protection for radio transmission operation. It may be necessary for the user to take measures, e.g., choose a different location for the device or realign it.

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Table 1
Guidelines and manufacturer's declaration – electromagnetic immunity

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment specified below. The user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ should make sure that it is used in such an environment.

Immunity tests	EN/IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) acc. to IEC 61000-4-2	± 8 kV contact discharge ± 15 kV air discharge	± 8 kV contact discharge ± 15 kV air discharge	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Rapid transient electrical Interferences/bursts acc. to IEC 61000-4-4	± 2 kV for power lines ± 1 kV for input and output lines 100 kHz repetition	± 2 kV/1kV ± 1 kV for input and output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surges acc. to IEC 61000-4-5	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	± 1 kV line(s) to line(s) ± 2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and fluctuations of the supply voltage acc. to IEC 61000-4-11	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	<u>Voltage Dip:</u> Dip to 0% for 1 cycles @ 0° phase angle Dip to 70% for 25/30 cycles @ 0° phase angle Dropout to 0% for 0.5 cycles @ 0°, 45°, 90°, 135°, 180°, 225°, 270° & 315° phase angles <u>Voltage interruption:</u> 100% for 250/300 cycles	Mains power quality should be that of a typical commercial or hospital environment. If the user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ requires continued operation during power supply interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Magnetic field at the supply frequency (50/60 Hz) acc. to IEC 61000-4-8	30 A/m at 50Hz/60Hz	30 A/m at 50Hz/60Hz	If image distortion occurs, it may be necessary to position the video mediastinoscope 10973 HD together with the IMAGE1 S™ further from sources of power frequency magnetic fields or to install magnetic shielding. The power frequency magnetic field should be measured in the intended installation location to assure that it is sufficiently low.
Radiated, radio-frequency, electromagnetic field immunity test acc. to IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz *Refer table 2 for wireless Proximity RF field test levels	3V/m 80 MHz to 2.7 GHz	
Immunity to conducted disturbances, induced by radio-frequency fields acc. to IEC 61000-4-6	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	3Vrms on 150 kHz to 80MHz 1kHz 80% AM modulation 6Vrms in ISM bands	

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Table 3
**Guidance and manufacturer's declaration – electromagnetic immunity –
for equipment and systems that are not life-supporting**

The Video-Mediastinoscope 10973 HD together with the IMAGE1 S™ is designed for operation in an electromagnetic environment as defined below. The user of the device should make sure that it is used in such an environment.

Interference immunity tests	EN/IEC 60601 test level	Conformity level	Electromagnetic environment – guidelines
Conducted RF disturbances acc. to IEC 61000-4-6	3 V _{rms} 150 kHz to 80 MHz	3 V _{rms}	Portable and mobile RF communications equipment should be used no closer to any part of the Video-Mediastinoscope 10973 HD together with the IMAGE1 S™, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance: $d = 1.2 \sqrt{P}$ where P is the nominal power of the transmitter in watts [W] according to the information provided by the manufacturer of the transmitter and d is the recommended separation distance in meters [m]. Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^a should be less than the $d = 1.2 \sqrt{P} \quad 80 \text{ MHz to } < 800 \text{ MHz}$ $d = 2.3 \sqrt{P} \quad 800 \text{ MHz to } 2.5 \text{ GHz}$ compliance level in each frequency range ^b . Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF disturbances acc. to IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

Note 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Note 2: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of buildings, objects and persons.

^a) Field strengths from fixed transmitters, such as base stations for radio [cellular/cordless] telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to stationary transmitters, an electromagnetic site survey should be considered. If the measured field strength at the location where the video mediastinoscope 10973 HD together with the IMAGE1 S™ is used exceeds the above compliance levels, the video mediastinoscope 10973 HD together with the IMAGE1 S™ should be monitored to ensure proper function. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the device.

^b) Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.

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Table 4

Guidelines and manufacturer's declaration – electromagnetic emissions

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in the electromagnetic environment specified below. The customer or user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ should ensure that it is used in such an environment.

RF emissions as per CISPR 11	Group 1	The video mediastinoscope 10973 HD together with the IMAGE1 S™ uses RF energy only for its internal function. The customer or user of the flexible video mediastinoscope 10973 HD together with the IMAGE1 S™ should make sure that it is used in such an environment. Therefore, its RF emissions are very low and are not likely to cause interference in nearby electronic equipment.
RF Emissions as per CISPR 11	Class A	The video mediastinoscope 10973 HD together with the IMAGE1 S™ is suitable for use in all establishments other than domestic and those directly connected to the public low voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions as per IEC 61000-3-2	Class A	
Voltage fluctuations/flicker Emissions as per IEC 61000-3-3	Complies	

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Table 5
Recommended separation distances between portable and mobile RF telecommunication devices and the Video-Mediastinoscope 10973 HD together with the IMAGE1 S™

The video mediastinoscope 10973 HD together with the IMAGE1 S™ is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or user of the video mediastinoscope 10973 HD together with the IMAGE1 S™ can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the video mediastinoscope 10973 HD together with the IMAGE1 S™ as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter [W]	Separation distance d [m] as a function of the transmitter frequency		
	150 kHz to 80 MHz $d = 1,2 \sqrt{P}$	80 MHz to 800 MHz $d = 1,2 \sqrt{P}$	800 MHz to 2.5 GHz $d = 2,3 \sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at a maximum output power not listed above, the recommended separation distance D in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for higher frequency range applies.

NOTE: These guidelines may not apply in all situations. The propagation of electromagnetic quantities is affected by absorptions and reflections of structures, objects and people.

HF surgical instrument compatibility

This equipment has been evaluated for compatibility with high-frequency surgical equipment.

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KARL STORZ Endoscopy Canada Ltd.
7171 Millcreek Drive, Mississauga, Ontario
L5N 3R3 Canada
Phone: +1 905 816-4500, Fax: +1 905 816-4599
Toll free phone: 1-800-268-4880 (Canada only)
Toll free fax: 1-800-482-4198 (Canada only)
E-Mail: info-canada@karlstorz.com

KARL STORZ Endoscopy-America, Inc.
2151 East Grand Avenue
El Segundo, CA 90245-5017, USA
Phone: +1 424 218-8100, Fax: +1 424 218-8525
Toll free phone: 800 421-0837 (USA only)
Toll free fax: 800 321-1304 (USA only)
E-Mail: communications@kseas.com

KARL STORZ Veterinary Endoscopy-America, Inc.
1 South Los Carneros Road
Goleta, CA 93117, USA
Phone: +1 805 968-7776, Fax: +1 805 685-2588
E-Mail: info@karlstorzvet.com

KARL STORZ Endoscopia Latino-America, Inc.
815 N. W. 57th Avenue, Suite 480
Miami, FL 33126-2042, USA
Phone: +1 305 262-8980, Fax: +1 305 262-8986
E-Mail: info@ksela.com

KARL STORZ Endoscopia México S.A. de C.V.
Edificio Atlantic, Oficina 3G
Calle D e/ 1ra y 3ra
10400 Vedado, Havana, Cuba
Phone: +537 836 95 06, Fax: +537 836 97 76
E-Mail: kstorz.cuba@gmail.com

KARL STORZ Endoscopia México S.A. de C.V.
Av. Ejercito Nacional No. 453 Piso 2,
Colonia Granada, Alcala Miguel Hidalgo
C.P. 11520 Ciudad de México
Phone: +52 (55) 1101 1520
E-Mail: mx-info@karlstorz.com

KARL STORZ Marketing América Do Sul Ltda.
Rua Joaquim Floriano, nº. 413, 20º andar – Itaim Bibi,
CEP-04534-011 São Paulo, Brasil
Phone: +55 11 3526-4600, Fax: +55 11 3526-4680
E-Mail: br-info@karlstorz.com

KARL STORZ Endoscopia Argentina S.A.
Zufriategui 627 6º Piso, B1638 CAA - Vicente Lopez
Provincia de Buenos Aires, Argentina
Phone: +54 11 4718 0919, Fax: +54 11 4718 2773
E-Mail: info@karlstorz.com.ar

KARL STORZ Endoskopi Norge AS
Stamveien 1
1483 Hagan, Norway
Phone: +47 6380 5600, Fax: +47 6380 5601
E-Mail: post@karlstorz.no

KARL STORZ Endoskop Sverige AB
Storsåtragränd 14
127 39 Skärholmen, Sweden
Phone: +46 8 505 648 00
E-Mail: kundservice@karlstorz.se

KARL STORZ Endoscopy Suomi Oy
Taivaltie 5
01610 Vantaa, Finland
Phone: +358 (0)96824774, Fax: +358 (0)968247755
E-Mail: asiakaspalvelu@karlstorz.fi

KARL STORZ SE & Co. KG
Representation Office
Kestuo st. 59 / Lenktoji st. 27
08124 Vilnius, Lithuania
Phone: +370 5 272 0448
Mobile: +370 685 67 000
E-Mail: info-lt-lv@karlstorz.com

KARL STORZ Endoskopi Danmark A/S
Skovlytoften 33
2840 Holte, Denmark
Phone: +45 45162600, Fax: +45 45162609
E-Mail: marketing@karlstorz.dk

KARL STORZ Endoscopy (UK) Ltd.
415 Perth Avenue, Slough
Berkshire, SL1 4TQ, United Kingdom
Phone: +44 1753 503500, Fax: +44 1753 578124
E-Mail: info-uk@karlstorz.com

KARL STORZ Endoscopia Nederland B.V.
Displayweg 2
3821 BT Amersfoort, Netherlands
Phone: +31 (0)33 4545890
E-Mail: info-nl@karlstorz.com

KARL STORZ Endoscopy Belgium N.V.
Phone: +31 (0)33 4545890
E-Mail: info-be@karlstorz.com

KARL STORZ Endoscopia France S.A. S.
12, rue Georges Guynemer, Quartier de l'Europe
78280 Guyancourt, France
Phone: +33 1 30484200, Fax: +33 1 30484201
E-Mail: marketing-fr@karlstorz.com

KARL STORZ Endoskop Austria GmbH
Landstraßer Hauptstr. 148/1/G1
1030 Wien, Austria
Phone: +43 1 71 56 0470, Fax: +43 1 71 56 0479
E-Mail: storz-austria@karlstorz.at

KARL STORZ Endoscopia Ibérica S.A.
Parque Empresarial San Fernando
Edificio Munich – Planta Baja
28830 Madrid, Spain
Phone: +34 91 6771051, Fax: +34 91 6772981
E-Mail: info-es@karlstorz.com

KARL STORZ Endoscopia Italia S. r. l.
Via dell'Artigianato, 3
37135 Verona, Italy
Phone: +39 045 8222000, Fax: +39 045 8222001
E-Mail: info-ita@karlstorz.com

KARL STORZ Croatia d.o.o.
Capraška 6
10000 Zagreb, Croatia
Phone: +385 1 6406 070, Fax: +385 1 6406 077
E-Mail: info@karlstorz.hr

KARL STORZ Endoskopija d.o.o.
Cesta v Gorice 34b
1000 Ljubljana, Slovenia
Phone: +386 1 620 5880, Fax: +386 1 620 5882
E-Mail: pisarna@karlstorz.si

KARL STORZ Endoskop Polska Marketing Sp. z o.o.
ul. Bojkowska 47
44-100 Gliwice, Poland
Phone: +48 32 706 13 00, Fax: +48 32 706 13 07
E-Mail: info-pl@karlstorz.com

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TOO KARL STORZ Endoscopy Kasakhstan
Saryarka, 6, BC "Arman", off. 910
010000 Astana, Republic of Kazakhstan
Phone: +7 7172 552-549, 552-788, Fax: -444
E-Mail: info@karlstorz.kz

KARL STORZ ENDOSKOPE
East Mediterranean & Gulf (branch)
Building West Side 7A – Unit 7WA – 3008
Dubai Airport Free Zone, P.O. Box 54983
Dubai – United Arab Emirates
Phone: +971 (0)4 2958887, Fax: +971 (0)4 3205282
Service Hotline: +971 (0)4 3415882
E-Mail: info-gne@karlstorz-emg.com

KARL STORZ Endoscopy India Private Limited
11th Floor, Dr. Gopal Das Bhawan
28, Barakhamba Road
New Delhi 110001, India
Phone: +91 11 4374 3000, Fax: +91 11 4374 3010
E-Mail: corporate@karlstorz.in

KARL STORZ SE & Co. KG
Interchange 21 Tower, Level 33
399 Sukhumvit Road,
North Klongtoey, Wattana,
10110 Bangkok, Thailand
Phone: +84 28 3823 8000
Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ SE & Co. KG
Resident Representative Office
14th Floor, M Plaza Saigon
39 Le Duan, District 1
Ho Chi Minh City, Vietnam
Phone: +84 28 3823 8000, Fax: +84 28 3823 8039
E-Mail: infovietnam@karlstorz.com

KARL STORZ Endoscopy China Ltd.
Room 2503-05, 25F AXA Tower, Landmark East,
No. 100 How Ming Street, Kwun Tong, Kowloon,
Hong Kong, People's Republic of China
Phone: +852 28 65 2411, Fax: +852 28 65 4114
E-Mail: inquiry@karlstorz.com.hk

KARL STORZ Endoscopy (Shanghai) Ltd., Beijing Branch
Room 1805-1807, Building B, 18F Beijing IFC
No. 8, Jianguomenwai Street, Chaoyang District,
100022, Beijing, People's Republic of China
Phone: +86 10 5638188, Fax: +86 10 5638199
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shanghai Branch
Room 701A Building 5 & Room 501 Building 7,
No. 3000 Longdong Avenue, Pilot Free Trade Zone,
201203, Shanghai, People's Republic of China
Phone: +86 21 60339888, Fax: +86 21 60339808
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Chengdu Branch
Room 803-805, 8F Jin Jiang International Building
No. 1 West Linjiang Road, Wuhou District,
6100414, Chengdu, People's Republic of China
Phone: +86 28 86587977, Fax: +86 28 86587975
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Shenyang Branch
Room 2001-2005, 20F N-MEDIA International Center,
No. 167 Youth Avenue, Shenhe District,
110014, Shenyang, People's Republic of China
Phone: +86 24 23181118, Fax: +86 24 23181119
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy (Shanghai) Ltd., Guangzhou Branch
Room 02B & 03 & 04A, 35F Teem Tower,
No. 208 Tianhe Road, Tianhe District,
510620, Guangzhou, People's Republic of China
Phone: +86 20 87321281, Fax: +86 20 87321286
E-Mail: info@karlstorz.com.cn

KARL STORZ Endoscopy Asia Marketing Pte Ltd.
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Singapore Sales Pte Ltd
No. 8 Commonwealth Lane #03-02
Singapore 149555, Singapore
Phone: +65 69229150, Fax: +65 69229155
E-Mail: infoasia@karlstorz.com

KARL STORZ Endoscopy Korea Co. Ltd.
9F Hyowon-Building
97 Jungdae-ro, Songpa-gu
05719 Seoul, Korea
Phone: +82-70-4350-7474, Fax: +82-70-8277-3299
E-Mail: infokorea@karlstorz.com

KARL STORZ Endoscopy Taiwan Ltd.
12F, No. 192, Sec. 2, Chung Hsin Rd.,
Sindian District, New Taipei City, Taiwan
Phone: +886 933 014 160, Fax: +886 2 8672 6399
E-Mail: info-tw@karlstorz.com

KARL STORZ SE & Co. KG
Representative Office Philippines
1901 Picadilly Star Bldg., 4th Avenue, BGC
Taguig City 1636, Philippines
Phone: +63 2 317 45 00, Fax: +63 2 317 45 11
E-Mail: philippines@karlstorz.com

KARL STORZ Endoscopy Japan K. K.
Stage Bldg. 8F, 2-7-2 Fujimi
Chiyoda-ku, Tokyo 102-0071, Japan
Phone: +81 3 6380-8622, Fax: +81 3 6380-8633
E-Mail: info@karlstorz.co.jp

KARL STORZ Endoscopy Australia Pty. Ltd.
68 Waterloo Road, Macquarie Park NSW 2113
P O Box 50 Lane Cove NSW 1595, Australia
Phone: +61 (0)2 9490 6700
Toll free: 1800 996 562 (Australia only)
Fax: +61 (0)2 9420 0695
E-Mail: karlstorz@karlstorz.com.au
www.karlstorz.com

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Risks of infection

The medical device is not sterile when delivered. The use of non-sterile medical devices poses a risk of infection for patients, users and third parties.

- Inspect medical devices for visible contamination. Reprocess contaminated medical devices.
- Reprocess medical devices before and after each use. When doing so, use validated procedures.

Improper use

Incorrect application of medical devices poses a risk of injury for patients.

- Users of medical devices must have an appropriate medical qualification and be acquainted with the application.

Risks due to damaged parts

Incorrectly assembled and damaged medical devices can lead to injuries to the patient.

- Check the following properties of medical devices before and after every use to ensure that they: are free from damage, are in full working order, have no unintentional rough surfaces, have no sharp corners, burred edges, have no projecting parts, are complete.
- Do not leave missing or broken-off components inside the patient.

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Maintenance and installation errors endanger the safety of physicians, medical support staff and patients. The medical devices doesn't work reliably or works with only with limited output.

- Only staff authorized by KARL STORZ may perform the following work: Assembly, Expansion, Readjustments, Modifications, Repairs.

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Primary processing at the point of use

Processing of devices should begin within two hours of their use to ensure the effectiveness of the processing processes specified in the processing instructions.

Wipe the surfaces of the device with a compress or a disposable cloth to remove coarse dirt, corrosive solutions and medicines.

Rinse the surfaces of the device with cold water.

Immediately after use / pre-treatment place the dry device in the shipping container.

Securely packaged device should be delivered to the place of processing.

Cleaning and disinfection

The following methods are validated and permitted for cleaning and disinfection of the device:

- Cleaning surfaces with a brush
- Mechanical cleaning: thermal disinfection

Processing with decontamination in a mechanized way

The product is not intended for ultrasonic processing.

Protective cap application

Place a protective cap on the plug to protect the plug contacts when cleaning and disinfecting.

Cleaning surfaces with a brush

Clean the surfaces of the device with a brush under running cold water.

Brush surfaces until all visible dirt is removed.

Rinse surfaces under cold running water.

Mechanical cleaning / thermal disinfection

A washer-disinfector for heat-resistant products should be used for this device. The washer-disinfector must comply with the requirements of DIN EN ISO 15883-1.

The A_0 value of the disinfection process must be observed.

When selecting a suitable drawer, the equipment manufacturer should be consulted to ensure that the medical device is flushed or flushed.

Place the device in the drawer and in the instrument holder if necessary.

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After disinfection the device should be completely dried using compressed air or an air-filled syringe.

Visual check

Remove the protective cap from the plug.

Check the product for the following:

- visible pollution;
- damage and corrosion;
- completeness;
- dryness.

Items with visible soiling should be subjected to a complete cleaning and disinfection process again.

Damaged or corroded medical devices must be discarded.

Incomplete device must be rejected or missing parts installed in them.

If necessary, dry the device by hand.

Function test

If the product does not meet one of the criteria below, or is visibly damaged, refer to the Maintenance, Repair and Disposal section in the instruction manual.

The following tests should be carried out to identify limitations in the functionality of the device:

- Check the surfaces of the medical device for mechanical damage or changes.
- Check the legibility of the labels.
- Check the product for mechanical damage.
- Check the correct fit of the assembled components and, if necessary, the cleaning connector.
- Check if the image is being transferred.
- Check if the light is being transmitted and that there is no dirt or damage on the light entry or exit surfaces.

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e-mail: Info-RU@karlstorz.com

Limitation of liability

KARL STORZ shall be liable for the safe operation, operational reliability, and performance of this equipment only subject to the following conditions:

- all assembly operations, system expansions, readjustments, modifications, or repairs have been performed by a person or persons duly authorized by KARL STORZ;
- all electrical installations at the location in which the device is connected and operated meet the applicable laws and standards; and
- the device has been used in accordance with its operating instructions at all times.

Manufacturer's warranty

The guarantees provided can be found in the Standard Conditions of Business of KARL STORZ.

The medical device must always be sent to your local subsidiary (see "Subsidiaries" section), even during the warranty period.

Opening the equipment or performance of any repairs or modifications to the equipment by unauthorized persons shall relieve us of any liability for its performance. Any such repair or modification performed during the warranty period shall void all warranty.

Warranty period: The standard warranty period for the device is 12 months from the date of commissioning of the device, but no more than 15 months from the date of shipment.

The procedure for submitting complaints and responses to them is governed by civil law. A complaint may be brought only on such issues that were not the subject of acceptance of the goods made in accordance with the terms of the contract.

For all issues related to product maintenance, please contact the Authorized Representative of the manufacturer.

In case of complaint, contact authorized representative of the manufacturer.

Authorized representative of the manufacturer:

Limited Liability Company "KARL STORZ - Endoscopy WOSTOK" (LLC "KSHTÉ VOSTOK")

115114 Moscow, Russia, Derbenyevskaya nab. 7, p. 4

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