«APPROVE»

Becton Dickinson and Company,
7 Loveton Circle, Sparks, Maryland 21152, USA

Date: November 28, 2016

Signature Carol C. Leford

position International Regulatory
Affairs Manager

Seal

Jennifer A. Creamer
11/28/2016

Jennifer A. Creamer
Notary Public, Maryland
Baltimore County
My Commission Expires.
October 19, 2019

**Operational Documentation
(User's Manual)
for Medical Device:**

BD Phoenix Automated Microbiology System, variations: BD Phoenix 100, BD
Phoenix M50, with Accessories

Manufacturer:

Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152, USA

BD Phoenix™

Automated Microbiology System User's Manual

2011/11
Document Number: L-003342(S)
 448021
Not for Sale in USA

Becton, Dickinson and Company
7 Loveton Circle
Sparks, Maryland 21152 USA
Tel: 800.638.8663

Benex Limited
Rineanna House
Shannon Free Zone
Shannon, County Clare, Ireland
Tel: 353.61.47.29.20
Fax: 353.61.47.29.07

Change History

Revision	Date	Pages	Reason
A	5/00	All	New product release
B	6/00	Preface (i, ii, vii); App. E (all); App. F (F-7 – F-13)	Update taxa
C	8/00	All	Update with clinical F/T data
D	10/00	i, ii, 7-6	Add CE marks and error subcode
E	2/01	i, ii, 8-7, 8-8, F-3, F-4	Add new drugs to Section 8, App. F
F	3/01	All	Add info. for software vers. 2
G	10/01	Preface (all); Sect. 2 (all); Sect. 3 (all); Sect. 4 (all); Sect. 5 (all); Sect. 6 (all); Sect. 8 (all); Glossary (all); Index (all); Reader Comment Card	Add info. for software vers. 3
H	1/02	i, ii, 8-5 – 8-8, F-5 – F-6	Correction to AM/S. aureus (App. F); add GMF and STF (Sect. 8.1.3)
I	3/02	i, ii, Sect. 2 (5, 16, 29, 31); 3-4; Sect. 4 (7, 8, 32); Sect. 5 (4, 5, 14); 6-2; 7-7; App. E (2, 6, 9); App. F (9 – 12)	Changes for software vers. 3.22
J	2003/06	All	Changes for software vers. 3.32
K	2003/11	Preface, Sect. 1, 2, 4, 5, 6, 7, 8 App. F, J, Index	Changes for software vers. 4.01
L	2004/10	All	Changes for software vers. 4.05
M	2005/03	Preface, Sect. 1 (10-12), Sect. 2, Sect. 3 (8), Sect. 4, Sect. 5, Sect. 6 (9, 16-21), Sect. 7 (5), Sect. 8, Glossary, App. E, App. F (12-22), App. J (1, 10, 12), Index	Changes for software vers. 5.02
N	2006/08	Preface (i-ii), Sect. 1 (1-8), Sect. 2 (11-12, 23-36), Sect. 4 (1-12, 23-24), Sect. 5 (15-18, 43-44), Glossary, App. E, App. F (9-24), App. H	Minor corrections; add info. on new claims
O	2008/01	All	Changes for software vers. 5.55
P	2008/09	Preface, Sect. 1 (2-4, 8), Sect. 2 (11, 15, 17, 20, 22), Sect. 4, Sect. 5 (39, 41, 43, 47-48, 53, 55, 60), Sect. 6 (6, 9, 10, 42 - 46), G-5, App. E, App. F, App. H, J-12, RCC	Changes for software vers. 5.66/ 4.61
Q	2010/09	Preface (i-ii), Sect. 1, Sect. 2 (8, 24, 26, 33, 34), Sect. 4 (4-8, 21), Sect. 5 (4, 22, 28, 47, 48, 52), Sect. 7 (6), Sect. 8, Glossary (2), App. E, App. F, App. H	Changes for PUD 4.91; correct typos and AR address
R	2011/05	Preface, Sect. 1 (1-3), Sect. 2 (12-13), Sect. 3 (10), Sect. 4 (all), Sect. 5 (all), Sect. 6 (9, 34-48), Sect. 7 (2, 18 - 22), Sect. 8, App. E, App. F (5, 13--28), App. H, Index	Changes for software vers. 6.01/ 5.11
S (19)	2011/11	i, ii, Sect. 8	To include the complete product name, update the performance for 2 gram-negative drugs and correct typos

BBL, BD, BDXpert, BD logo, Crystal, CrystalSpec, EpiCenter, Phoenix, and PhoenixSpec, are trademarks of Becton, Dickinson and Company. ATCC is a trademark of American Type Culture Collection. CHROMagar is a trademark of Dr. A. Rambach.

© Copyright Becton, Dickinson and Company, 2000 – 2011. All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written permission of BD, 7 Loveton Circle, Sparks, Maryland, 21152, United States of America.

Contents

1 – Introduction

1.1	Intended Use	1–1
1.2	Summary and Explanation of the Test	1–1
1.3	Principles of the Procedure	1–2
1.4	Precautions	1–6
1.5	System Overview	1–6
1.5.1	Instrument Overview	1–6
1.5.1.1	Instrument Optical System and Drive Mechanism	1–6
1.5.1.2	Carousel Assembly	1–6
1.5.1.3	Incubation System	1–7
1.5.1.4	Panel Status and Internal Barcode Scanner Assembly	1–7
1.5.1.5	Front-Panel and External Barcode Scanners	1–8
1.5.1.6	Control Electronics Assembly	1–8
1.5.2	Software and Operation Overview	1–8
1.5.3	Panel Overview	1–9
1.5.4	Testing Overview	1–10
1.5.5	Normalizers	1–10
1.6	Manual Structure	1–11
1.7	Use of this Manual	1–12
1.8	Conventions	1–12
1.8.1	Keys	1–12
1.8.2	Symbols Used on the Equipment	1–13
1.8.3	Notes, Cautions, and Warnings	1–14
1.9	Summary of Cautions and Warnings	1–14

2 – Installation

2.1	General	2–1
2.1.1	Installation	2–2
2.2	Instrument Specifications	2–2
2.3	Software Setup	2–4
2.3.1	Instrument Configuration	2–7
2.3.1.1	General Instrument Configuration	2–7
2.3.1.2	More Configuration Options	2–11
2.3.1.3	Password Administration	2–13

Contents

- 2.3.2 Communications Configuration 2-14
 - 2.3.2.1 Organism Configuration 2-19
 - 2.3.2.2 Antimicrobial Configuration 2-22
- 2.3.3 Custom Interpretation Rule Set 2-24
- 2.3.4 BDxpert Rule Configuration 2-32
- 2.3.5 Rapid Reporting Configuration 2-34
- 2.3.6 Panel Lot Definition 2-38
- 2.4 System Startup 2-43
- 2.5 Software Installation 2-43
- 3 – Controls and Indicators
 - 3.1 General 3-1
 - 3.2 Power Switch 3-3
 - 3.3 Keypad and LCD Display 3-3
 - 3.3.1 DOWN /UP ARROW Keys 3-3
 - 3.3.2 TAB Key 3-4
 - 3.3.3 Soft Keys 3-4
 - 3.3.4 UNLOAD PANELS Key 3-5
 - 3.3.5 LOAD PANELS Key 3-5
 - 3.3.6 SILENCE ALARM Key 3-5
 - 3.3.7 System Alert/Door Unlocked Indicator 3-5
 - 3.3.8 Display Area 3-5
 - 3.4 Barcode Scanners 3-6
 - 3.5 Floppy Disk Drive 3-6
 - 3.5.1 Floppy Disk Indicator 3-6
 - 3.5.2 Floppy Disk Eject Button 3-6
 - 3.6 Station Indicators 3-7
 - 3.7 Audible Tones and Alarms 3-8
 - 3.8 Keyboard 3-9
 - 3.9 Printer 3-9

Contents

4 – Operation

4.1	General.....	4-1
4.2	Daily System Maintenance.....	4-2
4.3	Storage and Handling.....	4-2
4.4	Preparing Panels.....	4-2
4.4.1	Quality Control	4-15
4.5	Panel Login.....	4-19
4.5.1	Using the Instrument Interface.....	4-19
4.5.2	Logging in Panels	4-19
4.5.3	Associating and Removing LIS Orders.....	4-21
4.6	Inserting Panels in the Instrument.....	4-24
4.7	Panel Testing	4-26
4.8	Responding to Needs Attention Panels.....	4-26
4.9	Automatic Association of Panels	4-30
4.9.1	Overview	4-30
4.9.2	Panel Types.....	4-31
4.9.3	When Auto Association Is Not Performed.....	4-31
4.9.4	Organism ID Conflict Checking	4-31
4.9.5	Auto Association and Related Actions	4-31
4.9.6	Typical Auto Association Examples	4-32
4.10	Obtaining Results	4-34
4.10.1	Adding/Modifying ID Results	4-38
4.10.2	Modifying AST Results	4-39
4.10.3	BDXpert Triggered Rules Display.....	4-40
4.10.4	Resistance Marker Display	4-40
4.11	Panel Finalization.....	4-41
4.12	Responding to Alarms and Errors	4-45
4.13	Printing Reports.....	4-46
4.14	Unloading and Discarding Panels	4-48
4.15	Installing New Panel Formats	4-49
4.16	Phoenix Update Disk (PUD).....	4-50

Contents

4.17	LIS Operations	4-50
4.17.1	General	4-50
4.17.2	Important Concepts	4-52
4.17.3	Processing Panel Orders	4-52
4.17.4	Routine System Operation.....	4-52
 5 – Reference		
5.1	General	5-1
5.2	Software Tree.....	5-2
5.3	Main Status Screen	5-3
5.4	Panel Login	5-6
5.5	Panel Results	5-10
5.5.1	BDXpert Triggered Rules Display	5-20
5.5.2	Resistance Marker Display	5-23
5.5.3	Special Messages Display	5-24
5.5.4	Lab Report	5-27
5.6	Finalization	5-29
5.7	Configuration/Maintenance	5-31
5.8	Reports	5-32
5.8.1	Completed Lab Report	5-33
5.8.2	Accession Lab Report	5-36
5.8.3	Needs Attention List Report	5-38
5.8.4	Resident Panel Report	5-39
5.8.5	QC Lab Report.....	5-41
5.8.6	Cumulative QC Report.....	5-43
5.8.7	Daily Instrument Report	5-45
5.8.8	Interpretation Rule Set Report.....	5-47
5.8.9	BDXpert Rule Set Database Report.....	5-50
5.8.10	Organism ID Code List Report	5-52
5.8.11	Antimicrobial Code Report	5-53
5.8.12	Lab Report/QC Lab Report.....	5-54
5.8.13	Finalization Summary Report	5-54
5.8.14	Custom Breakpoint Difference Report	5-55

Contents

5.8.15	Current QC Panel Lot Report.....	5-55
5.8.16	Historical QC Panel Lot Report.....	5-55
5.8.17	Panel Lot Report.....	5-55
5.8.18	Panel Lot Database Report.....	5-57
5.9	Panel Inventory	5-58
5.10	System Alerts.....	5-61
5.11	Needs Attention	5-62
6 – Maintenance		
6.1	General.....	6-1
6.2	Routine Maintenance	6-2
6.2.1	Daily Maintenance	6-2
6.2.2	Weekly Maintenance	6-2
6.2.3	Periodic Maintenance	6-4
6.2.3.1	Air Filter Replacement.....	6-4
6.2.3.2	Cleaning the Barcode Scanner Window	6-7
6.2.3.3	Cleaning/Decontamination	6-8
6.2.4	Maintenance Menu	6-9
6.2.4.1	Save User Data to Disk.....	6-9
6.2.4.2	Restore User Data from Disk	6-11
6.2.4.3	Save LIS Codes to Disk	6-13
6.2.4.4	Restore LIS Codes from Disk.....	6-14
6.2.4.5	Save Panel Lot Definitions	6-16
6.2.4.6	Restore Panel Lot Definitions	6-18
6.2.4.7	Save Panel Configuration	6-20
6.2.4.8	Install Panel Configuration.....	6-22
6.2.4.9	Phoenix Update Disk	6-24
6.2.4.10	Modify Panel Usage	6-25
6.2.4.11	Write Light Source to Disk	6-27
6.2.4.12	Read Light Source from Diskette.....	6-29
6.2.4.13	View Syslog.....	6-31
6.2.4.14	Save Syslog to Network	6-34
6.2.4.15	Display UV and RGB Status	6-36
6.2.4.16	Force Front UV Bulb Adjustment.....	6-37
6.2.4.17	Force Rear UV Bulb Adjustment	6-39
6.2.4.18	Force RGB LED Adjustment	6-42
6.2.4.19	Front Panel Reader Setup.....	6-44
6.3	Module Replacement.....	6-45
6.3.1	General.....	6-45
6.3.2	Barcode Scanner Replacement.....	6-46
6.3.3	Thermometer Removal	6-47
6.3.4	Reuniting Separated Liquid in the Thermometer	6-47

Contents

7 – Troubleshooting

7.1	General	7-1
7.1.1	Instrument Service	7-1
7.2	Error/Alert Messages	7-2
7.3	Save System Data to Disk Function	7-27
7.4	Save Syslog to Disk Function	7-28

8 – Performance and Limitations

8.1	Performance Characteristics	8-1
8.1.1	Gram Negative Performance	8-1
8.1.2	Gram Positive Performance	8-4
8.1.3	Streptococci Performance (with the Phoenix SMIC/ID, SMIC Panels)	8-7
8.1.4	Identification of Yeast Species	8-9
8.2	Limitations of the Procedure	8-9

G – Glossary and Abbreviations

A – Limited Warranty

B – Replacement Parts

C – Software Update Log

C –

D – International Contacts

E – Codes

E.1	Organism Codes, Short Names, Long Names	E-1
-----	---	-----

Contents

F – Panel Information

F.1	List of Antimicrobial Agents in Phoenix Panels	F-1
F.2	List of Reagents and Principles Employed in the Phoenix System.....	F-8
F.2.1	Gram Negative	F-8
F.2.2	Gram Positive	F-9
F.2.3	Streptococci Panel	F-11
F.2.4	Yeast Panel	F-13
F.3	Taxa for ID/AST Determination	F-15
F.3.1	Gram Negative (0.5 McFarland).....	F-15
F.3.2	Gram Negative (0.25 McFarland).....	F-19
F.3.3	Gram Positive (0.5 McFarland).....	F-22
F.3.4	Gram Positive (0.25 McFarland).....	F-24
F.3.5	Streptococci.....	F-26
F.3.5	Yeast.....	F-27

H – Bibliography

J – Event Log Messages

Notes –

Notes – Notes

I – Index

2 – Reader Comment Card

Contents

Figures

Figure 1-1 – BD Phoenix 100 Instrument.....	1-7
Figure 1-2 – Panel Inoculation Station	1-9
Figure 1-3 – Panel Transportation Caddy	1-10
Figure 1-4 – Symbols Used on the Phoenix Instrument	1-13
Figure 2-1 – Main Status Screen	2-4
Figure 2-2 – Configuration Menu.....	2-6
Figure 2-3 – Instrument Configuration Display	2-9
Figure 2-4 – More Configuration Options Display	2-12
Figure 2-5 – Password Administration Display	2-14
Figure 2-6 – Communications Configuration Display	2-17
Figure 2-7a – Organism Configuration Display	2-20
Figure 2-7b – Organism Edit Display	2-21
Figure 2-8a – Antimicrobial Configuration Display.....	2-23
Figure 2-8b – Antimicrobial Edit Display	2-23
Figure 2-9a – Custom Interpretation Rule Set – Screen A.....	2-28
Figure 2-9b – Custom Interpretation Rule Set – Screen B	2-28
Figure 2-10 – Custom Breakpoint Difference Report	2-30
Figure 2-11 – BDxpert Rule Set Configuration	2-33
Figure 2-12 – Rapid Reporting Configuration Display	2-36
Figure 2-13 – Panel Lot Definition Display	2-39
Figure 2-14 – Sample Current QC Panel Lot Report.....	2-41
Figure 2-15 – Sample Historical QC Panel Lot Report	2-42
Figure 3-1 – Phoenix Instrument Layout	3-2
Figure 3-2 – On/Off Switch	3-3
Figure 3-3 – Keypad and LCD Display.....	3-4
Figure 3-4 – Floppy Disk Drive Controls and Indicators	3-6
Figure 3-5 – Station Status Indicators	3-7
Figure 4-1 – Panel and Closures	4-6
Figure 4-2 – QC Panel Login Display.....	4-17
Figure 4-3 – Panel Login – LIS Orders Display.....	4-23
Figure 4-4 – Inserting Panels	4-25
Figure 4-5 – Panel Results Display – ID Results	4-36
Figure 4-6 – Panel Results Display – AST Results	4-37
Figure 4-7 – Panel Results Display – ID Results Showing Triplets	4-39
Figure 4-8 – Resistance Marker Review Display	4-41
Figure 4-9 – Finalization Display	4-43
Figure 4-10 – Finalization Summary Report	4-44
Figure 4-11 – Reports Menu	4-47
Figure 5-1 – Main Status Screen (with optional LIS indicator).....	5-6
Figure 5-2 – Panel Login Display	5-7
Figure 5-3 – Panel Results List Display.....	5-13
Figure 5-4a – Panel Results Display – ID	5-15
Figure 5-4b – Panel Results Display – AST	5-17
Figure 5-5 – BDxpert Rules Triggered Display	5-22
Figure 5-6 – Resistance Marker Review Display	5-24
Figure 5-7 – Special Messages Display	5-26

Figures

Figure 5-8 – Sample Lab Report.....	5-28
Figure 5-9 – Finalization Display	5-31
Figure 5-10 – Sample Completed Lab Report	5-35
Figure 5-11 – Sample Accession Lab Report	5-37
Figure 5-12 – Sample Needs Attention List Report	5-39
Figure 5-13 – Sample Resident Panel Report	5-40
Figure 5-14 – Sample QC Lab Report.....	5-42
Figure 5-15 – Sample Cumulative QC Report.....	5-44
Figure 5-16 – Sample Daily Instrument Report	5-46
Figure 5-17 – Sample Interpretation Rule Set Report (sample data only).....	5-48
Figure 5-18 – Sample BDxpert Rule Set Database Report	5-51
Figure 5-19 – Sample Organism ID Code List Report	5-52
Figure 5-20 – Sample Antimicrobial Code Report	5-54
Figure 5-21 – Sample Panel Lot Report.....	5-56
Figure 5-22 – Sample Panel Lot Database Report	5-58
Figure 5-23 – Panel Inventory Display (primary)	5-59
Figure 5-24 – Panel Inventory Display (secondary).....	5-59
Figure 5-25 – System Alerts Display.....	5-61
Figure 5-26 – Needs Attention Display	5-65
Figure 6-1 – Main Status Screen	6-2
Figure 6-2 – Electronics Bay Filter Replacement.....	6-5
Figure 6-3 – Source Bay Filter Replacement	6-6
Figure 6-4 – Incubator Filter Replacement.....	6-7
Figure 6-5 – Barcode Scanner Replacement.....	6-46

Tables

2-1	Physical Dimensions	2-2
2-2	Environmental Requirements	2-3
2-3	Electrical Requirements	2-3
2-4	Optical Specifications	2-3
2-5	Date Formats	2-10
2-6	Rules Updates	2-25
3-1	Station Indicators	3-7
3-2	Audible Tones/Alarms	3-8
4-1	Recommended Media	4-3
4-2	Needs Attention Resolutions	4-29
5-1	MIC/SIR Values and Causes	5-19
5-2	Organism Group Codes	5-49
5-3	Needs Attention: Active Soft Keys	5-66

1 Introduction

1.1 Intended Use

The BD Phoenix™ Automated Microbiology System is intended for the rapid identification (ID) and antimicrobial susceptibility testing (AST) of clinically significant bacteria. The Phoenix system provides rapid results for most aerobic and facultative anaerobic Gram-positive bacteria as well as most aerobic and facultative anaerobic Gram-negative bacteria of human origin. The Phoenix system is also intended for the rapid identification of yeast and yeast-like organisms.

1.2 Summary and Explanation of the Test

Micromethods for the biochemical identification of microorganisms were reported as early as 1918¹. Several publications reported on the use of the reagent-impregnated paper discs and micro-tube methods for differentiating enteric bacteria²⁻⁹. The interest in miniaturized identification systems led to the introduction of several commercial systems in the late 1960s, and they provided advantages in requiring little storage space, extended shelf life, standardized quality control, and ease of use.

Many of the tests used in the Phoenix ID panels are modifications of the classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition to these, the Phoenix System utilizes chromogenic and fluorogenic substrates as well as single carbon source substrates in the identification of organisms.^{10, 11, 26, 27, 28}

The modern broth microdilution test used today has origins in the tube dilution test used in 1942 by Rammelkamp and Maxon to determine *in vitro* antimicrobial susceptibility testing of bacterial isolates from clinical specimens.¹² The broth dilution technique involves exposing bacteria to decreasing concentrations of antimicrobial agents in liquid media by serial two-fold dilution. The lowest concentration of an antimicrobial agent in which no visible growth occurs is defined as the minimal inhibitory concentration (MIC).

The introduction in 1956 of a microtitrator system, using calibrated precision spiral wire loops and droppers for making accurate dilutions rapidly allowed Marymont and Wentz to develop a serial dilution AST test.¹³ The microtitrator system was accurate and allowed the reduction in volumes of antimicrobial agents. The term microdilution appeared in 1970 to describe the MIC tests performed in volumes of 0.1 ml or less of antimicrobial solution.¹⁴

The Phoenix AST is a modified miniaturized version of the micro-broth doubling dilution technique. Susceptibility testing in the Phoenix System is performed through determination of bacterial growth in the presence of various concentrations of the antimicrobial agent tested with the aid of the AST indicator in continuously incubated and read micro-wells in the Phoenix panels.

1.3 Principles of the Procedure

A maximum of 100 identification and antimicrobial susceptibility tests can be performed in the Phoenix instrument at a time using Phoenix combination panels. A sealed and self-inoculating molded polystyrene tray, with 136 micro-wells containing dried reagents, serves as the Phoenix disposable. The combination panel includes an ID side with dried substrates for bacterial or yeast identification and an AST side with varying concentrations of antimicrobial agents, growth and fluorescent controls at appropriate well locations. The Phoenix system utilizes an optimized colorimetric redox indicator for AST, and a variety of colorimetric and fluorometric indicators for ID. The AST broth is cation-adjusted (e.g., Ca⁺⁺ and Mg⁺⁺) to optimize susceptibility testing performance.

The Phoenix Panel is comprised of a 51 well ID side and an 85 well AST side. The ID side contains 45 wells with dried biochemical substrates and 2 fluorescent control wells. The AST side potentially contains up to 84 wells with dried antimicrobial agents and 1 growth control well. Panels are available as ID only, AST only, or ID/AST combination. Unused wells are reserved for future use.

Phoenix panels are inoculated with a standardized inoculum. Organism suspensions must be prepared only with the BBL™ CrystalSpec™ nephelometer, the BD PhoenixSpec™ nephelometer, or the BD Phoenix AP™ instrument. Once inoculated, panels are placed into the instrument and continuously incubated at 35° C. The instrument tests panels every 20 minutes: on the hour, at 20 minutes past the hour, and again at 40 minutes past the hour for up to 16 hours if necessary. Phoenix panels are read only by the Phoenix instrument. Phoenix panels cannot be read manually.

Organism Identification

The ID portion of the Phoenix panel utilizes a series of conventional, chromogenic, and fluorogenic biochemical tests to determine the identification of the organism. Both growth-based and enzymatic substrates are employed to cover the different types of reactivity within the range of taxa. The tests are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Acid production is indicated by a change in phenol red indicator when an isolate is able to utilize a carbohydrate substrate. Chromogenic substrates produce a yellow color upon enzymatic hydrolysis of either P-nitrophenyl or P-nitroanilide compounds. Enzymatic hydrolysis of fluorogenic substrates results in the release of a fluorescent coumarin derivative. Organisms that utilize a specific carbon source reduce the resazurin-based indicator. In addition, there are other tests that detect the ability of an organism to hydrolyze, degrade, reduce, or otherwise utilize a substrate.

A complete list of taxa that comprises the Phoenix Database is provided in Appendix F, Section F.3. Reactions employed by various substrates and a brief explanation of the principles employed in the Phoenix Gram Negative, Gram Positive, *Streptococcus*, and Yeast ID reactions are described in Appendix F, Section F.2.

Antimicrobial Susceptibility Testing

The Phoenix AST method is a broth based microdilution test. The Phoenix system utilizes a redox indicator for the detection of organism growth in the presence of an antimicrobial agent.¹⁵ Continuous measurements of changes to the indicator as well as bacterial turbidity are used in the determination of bacterial growth. Each AST panel configuration contains several antimicrobial agents with a wide range of two-fold doubling dilution concentrations. Organism identification is used in the interpretation of the MIC values of each antimicrobial agent. A complete list of taxa for which the Phoenix panels can provide AST results is shown in Appendix F, Section F.3. The list of antimicrobial agents and concentrations available for susceptibility testing in the Phoenix system is provided in Appendix F, Section F.1.

The components required for testing include 1) Phoenix panels with panel closures, 2) Phoenix ID Broth, 3) Phoenix AST Broth, 4) Phoenix AST Indicator solution, 5) Phoenix AST-S Broth (for use with Phoenix Strep panels only), 6) Phoenix AST-S Indicator Solution (for use with Phoenix Strep panels only), 7) Phoenix Inoculation Station, 8) Phoenix Panel Caddy, 9) the BBL™ CrystalSpec™ nephelometer, the BD PhoenixSpec™ nephelometer, or the BD Phoenix AP™ instrument, 10) 25 µL pipettor and sterile tips, and 11) Miscellaneous lab supplies (listed under Section 4.4 – Materials Required but Not Provided).

The Phoenix disposable is for use with the Phoenix instrument for either organism identification (ID) or antimicrobial susceptibility testing (AST), or for the combination of both. The panel is inclined with the inoculation ports at the top for filling. Separate inocula are added manually to the ID and AST ports. The inocula flow down the panel in serpentine fashion, filling the panel wells as the liquid front progresses toward the pad. The pad absorbs excess inoculum. Closures are manually inserted in the fill ports. An air admittance port is located in the divider area of the panel lid to ensure adequate oxygen tension in the panel for the duration of the test.

Principles of Phoenix AST Tests for the Detection of Resistance Markers

The following sections outline the principles of the Phoenix AST System in the detection of resistance markers in gram-negative or gram-positive organisms, including 1) detection of ESBL production among species of *Enterobacteriaceae*; 2) detection of vancomycin resistance in *Enterococcus* species (VRE); 3) detection of high-level aminoglycoside resistance in *Enterococcus* and *Streptococcus* species (HLAR); 4) detection of methicillin-resistance in staphylococci (MRS); 5) detection of β-lactamase production in *Staphylococcus* species (BL); 6) detection of macrolide resistance in *Streptococcus* species (MLSb); 7) detection of mecA-mediated Resistance with *S. aureus* (mecA); 8) detection of Vancomycin Resistant *Staphylococcus aureus* (VRSA); 9) detection of Phoenix Inducible Macrolide Resistance (iMLSb); 10) detection of Phoenix high level Mupirocin resistance (HLMUPH).

Phoenix Extended Spectrum β-Lactamase (ESBL) Test¹⁶

The Phoenix ESBL test evolved from published data of known ESBL antibiogram patterns in the current literature.^{18,19,20,21} Selected strains of various species with known β-lactamase genotype/phenotypes in the family *Enterobacteriaceae*, including *Escherichia coli*, *Klebsiella* species (spp.), *Citrobacter* spp., *Enterobacter* spp., *Proteus* spp., and *Serratia* spp., were used to develop the Phoenix ESBL test. The Phoenix ESBL test is based on the principle of a differential response between the inhibitory effect of selected second or third generation cephalosporins in the presence or absence of a β-lactamase inhibitor, clavulanic acid. The principles of Phoenix ESBL test is similar to the CLSI ESBL broth microdilution confirmatory test²². The Phoenix ESBL test is applied to *E. coli*, *K. pneumoniae* and *K. oxytoca*. Additionally, at the customers discretion, it can be applied to other enteric species where production of ESBL has been reported in literature. When a test result of ESBL is positive, the categorical interpretation of all penicillins, cephalosporins (except cephamycins), and aztreonam on the same Phoenix panel will be changed to resistant.²³

Phoenix Vancomycin Resistant Enterococci (VRE) Test

The Phoenix VRE test is based on the SIR interpretation of vancomycin. The breakpoint selected in the instrument configuration is used for the categorical interpretation. The Phoenix VRE test was developed and optimized to match the CLSI standard broth microdilution test.^{22,23} Selection of a breakpoint other than CLSI may result in less than optimal performance due to differences in categorical interpretations. Only *Enterococcus faecalis* and *E. faecium* with acquired resistance (vanA or vanB) will be reported as positive.²²

Phoenix High-Level Aminoglycoside Resistance (HLAR) Tests

The Phoenix HLAR tests for *Enterococcus* are based on the growth response in a single well containing either a high-level concentration of gentamicin or streptomycin. These tests were developed and optimized against both the CLSI standard broth microdilution and the CLSI screening agar test.²²

The Phoenix HLAR tests for *Streptococcus* are based on the growth response in a single well containing gentamicin, kanamycin, or streptomycin. These tests were developed and optimized using the CLSI recommended standard broth microdilution.

Phoenix Methicillin-Resistance in Staphylococci (MRS) Test

The Phoenix MRS test is based on the SIR interpretation of oxacillin with *Staphylococcus* species. The breakpoint selected in the instrument configuration is used for the categorical interpretation. This test was developed and optimized against the CLSI standard broth microdilution^{22,23} and the CLSI oxacillin-salt agar screening test (BBL™ Oxacillin Screen Agar)²². Selection of a break-point other than CLSI may result in less than optimal performance due to differences in categorical breakpoints between standards. All *Staphylococcus* species will be evaluated with the Phoenix MRS test. When a test result of MRS is positive, the interpretations of all beta-lactam drugs on the same Phoenix panel will be changed to resistant.²²

Phoenix Gram-Positive β -lactamase (BL) Test¹⁶

The BL test available in the Phoenix AST System is a nitrocefin based β -lactamase test. The nitrocefin based test is a direct detection method located on the ID side of the Phoenix panel. The performance of this test was established against the results of testing with Cefinase™ Discs (BBL™, Cat. No. 31650) as the reference method. Currently, only *Staphylococcus* species will be evaluated with these tests. When the result of BL test is positive, the categorical interpretation of all penicillinase labile penicillins on the same Phoenix panels will be changed to resistant.²²

Phoenix Macrolide Resistance in Streptococci (MLSb) Test

The Phoenix Macrolide Resistance test is based on SIR interpretation of erythromycin and clindamycin. The breakpoint selected in the instrument configuration is used for the categorical interpretation. Erythromycin resistant and clindamycin resistant *Streptococcus* isolates will be reported as macrolide/lincosamide/streptogramin B (MLSb) phenotype.

Phoenix *mecA*-mediated Resistance Marker for *Staphylococcus aureus* (*mecA*)

The Phoenix *mecA* test is used to predict *mecA*-mediated resistance in *Staphylococcus aureus*. The principle is similar to the CLSI-recommended Disk Diffusion test, which uses a cefoxitin (FOX) disk to predict *mecA*-mediated resistance in *S. aureus*. The performance of the test was established against multiplex PCR methods²⁵ as well as the Disk Diffusion test. With the Phoenix *mecA* test, the *mecA*-specific FOX breakpoints used for detection of the resistance marker will be configured in the instrument. When the *mecA* resistance marker is detected, the interpretations for all beta-lactam drugs on the same Phoenix panel are changed to resistant 22, and the Phoenix *mecA* resistance marker is set.

Phoenix Vancomycin Resistant *Staphylococcus aureus* (VRSA) Test

The Phoenix VRSA detection is based on the SIR interpretation of vancomycin when testing *Staphylococcus aureus*. The breakpoint selected in the instrument configuration is used for the categorical interpretation. The Phoenix VRSA test was developed and optimized to match the CLSI standard broth microdilution test, and verified with known VRSA isolates. Selection of a breakpoint other than those found in CLSI M100-S15 may result in less than optimal performance due to differences in categorical interpretations. Only *Staphylococcus aureus* with true resistance (isolates containing resistance marker such as *vanA* gene) will be reported as VRSA. Strains of *S. aureus* with vancomycin intermediate results (GISA/VISA) will be identified and reported by separate BDxpert rules. The Phoenix Gram Positive AST panel detected vancomycin resistance in the VRSA *S. aureus* strains available at the time of comparative testing. The ability to detect resistance in other *S. aureus* strains is unknown due to the limited number of resistant strains available for comparative testing.

Phoenix Inducible Macrolide Resistance (iMLSb) Test in *Staphylococcus* species

The Phoenix Inducible Macrolide Resistance (iMLSb) Test is used to detect inducible macrolide-lincosamide-streptogramin B (MLSb) resistance in *Staphylococcus* species. MLSb resistance, usually encoded by *ermA* or *ermC* genes, may be either constitutive (always expressed) or inducible after exposure to a macrolide antibiotic (e.g. erythromycin, clarithromycin, etc.). The Phoenix Inducible Macrolide Resistance Test is based on the same principle as the CLSI-recommended Disk Approximation Test (D-Test) for the detection of inducible clindamycin resistance. When the Phoenix Inducible Macrolide Resistance Test result is positive, the categorical interpretation of clindamycin on the same Phoenix panel will be reported as resistant and accompanied by a separate BDxpert message. *Staphylococcus* isolates resistant to both erythromycin and clindamycin on initial testing will be reported as constitutive MLSb resistance to distinguish them from isolates that are resistant to macrolides alone by efflux mechanism.

Phoenix high level Mupirocin resistance (HLMUPH)

The Phoenix high level Mupirocin resistance test is based on the growth response in a single well containing 256 mcg/mL of the topical antibiotic Mupirocin that is frequently used for eradication of *Staphylococcus* colonization, particularly for nasal carriage of MRSA. The breakpoint selected in the instrument configuration is used for the categorical interpretation of Susceptible (no growth in 256 mcg/mL well) or Resistant (growth in 256 mcg/mL well). High level Mupirocin resistant *Staphylococcus* isolates will be reported as the (HLMUPH) phenotype resistance marker.

1.4 Precautions

For *in vitro* Diagnostic Use.

All patient specimens and microbial cultures are potentially infectious and should be treated with universal precautions. Please refer to CDC manual *Bio-safety in Microbiological and Biomedical Laboratories*, 4th Edition, 1999, as well as other recommended literature.

Panels once inoculated should be handled carefully until placed in the instrument.

1.5 System Overview

1.5.1 Instrument Overview

1.5.1.1 Instrument Optical System and Drive Mechanism

Within the main instrument incubation bay, the cylindrical carousel and its drive are the only moving parts. The carousel assembly is supported and permitted to rotate by three large V-groove bearing assemblies, which accept the carousel V-ring. The carousel is normally driven at either 1.0 RPM or 2.0 RPM. The drive speed, acceleration and deceleration are controlled by the central processor.

The carousel is divided vertically into four tiers, each of which functions as a relatively independent optical source and detection system for panels placed on that tier. Each tier has its own microcontroller to control data acquisition and transmission. The four microcontrollers communicate with the central processor over a serial communications link.

Visible illumination in the red, green and blue spectral regions for each tier is provided by a Light Emitting Diode (LED) source board. The LED source currents are programmable in order to compensate for signal loss at the panel extremes due to parallax and other factors. A source monitor system averages signal from two onboard photodiodes to monitor visible source output. A fluorescent UV (ultraviolet) lamp containing a selected phosphor provides UV illumination. The active length of the lamp spans all four tiers vertically, producing UV illumination constant to 10% over the useful lamp length. A second UV lamp is maintained as a spare.

1.5.1.2 Carousel Assembly

The Carousel Assembly is a cage-like structure comprised of aluminum rings and vertical ribs bolted together to form a right circular cylinder. The carousel holds 104 panel holders in four tiers. Each tier holds one normalizer panel and 25 sample panels.

During normal operation, the carousel is driven at either 1.0 RPM or 2.0 RPM, depending on the current operation. During panel location, the carousel can rotate at up to 10.0 RPM. The carousel rotates counterclockwise. A single rotation of the carousel is used to perform an inventory scan to identify panels within the instrument following a door closure. A complete test cycle requires seven minutes.


Figure 1-1 – BD Phoenix 100 Instrument

1.5.1.3 Incubation System

The Incubation System maintains the carousel, panel carriers, and panels at a constant nominal temperature of 35° C. The system is a recirculating forced-air convection design, and has a single filter to remove dust. The system consists of a squirrel cage blower powered by a brushless DC motor, a coil-style heater section with automatic over-temperature shut-off, thermoformed inlet and return ductwork, and a user-replaceable polyester fiber air filter.

1.5.1.4 Panel Status and Internal Barcode Scanner Assembly

Panel status indication and panel identification is accomplished on the same tower in the instrument, mounted directly behind the carousel in the door area. Panel Status is indicated by tri-color LEDs mounted on a panel behind the carousel that shine through light pipes in the panel carriers. Five columns of the carousel are exposed when the instrument door is open; there are five LED boards with four indicator LEDs on each.

Four barcode scanners capable of reading panel barcodes are located on the same tower as the indicators. Each scanner reads panel barcodes on its respective tier of the carousel. Panel barcodes are affixed at the top of the panel on the base (bottom) side of the panel, facing inward on the carousel when placed in the instrument.

1.5.1.5 Front-Panel and External Barcode Scanners

External barcode scanners are located on the front panel and the right side of the instrument. These scanners can be used to read barcoded accession numbers you have placed on the panels, as well as the panel's own sequence number barcodes. The accession barcodes can be used to link specimen identification information to specific panels in the instrument.

1.5.1.6 Control Electronics Assembly

The control electronics assembly consists of the boards: I/O (Input/Output) Interface, CPU and memory, Ethernet, and LCD (Liquid Crystal Display) Driver. Other control electronics include user interface components such as the front panel barcode scanner, the floppy disk drive, the LCD display, the custom keypad, the keyboard, and the speaker for audible alarm notification.

The control electronics are primarily responsible for the following functions:

- scheduling and monitoring panel testing
- acquiring panel test data
- illuminating the panel (station) status indicators
- processing data from the external and internal barcode scanners
- interpreting data from organism identification panels to produce an organism identification
- interpreting data from antimicrobial susceptibility panels to produce either MIC results or, via guidelines, SIR results
- sending information to the LCD display
- processing commands from the keypad and keyboard
- controlling door switch status and door interlock
- controlling cabinet incubation
- communications with external devices (such as an EpiCenter™ workstation)
- generating alarm signals
- performing the Built-In-Test (BIT)

1.5.2 Software and Operation Overview

The Liquid Crystal Display (LCD) presents all the information needed to monitor instrument status, to enter and remove panels, set up the instrument, print reports, and perform routine instrument maintenance. The information is presented in the form of screens and icons that graphically represent the information (e.g., a clock indicates the current time). The top region of the Main Status display presents instrument status information that is updated every few seconds. The middle region of the display initially presents station status information. As you access different screens and functions, they appear in place of the top and middle status windows. At the bottom of the displays, a series of icons shows the current software (soft) key definitions. Each of the displays are discussed in detail in **Section 5 – Reference**.

Most of the operations you perform at the instrument are initiated by pressing soft keys. Soft key definitions change as you select their related functions and access different displays and operations.

In addition to the soft keys, six hard keys, themselves marked with screened icons, perform fixed functions regardless of the current display or operation. The LOAD PANELS key is used to prepare the instrument for loading panels for testing; the UNLOAD PANELS key is used to remove panels from the instrument; the SILENCE ALARM key quiets any audible alerts; the UP and DOWN ARROW keys are used to scroll through lists or to increase or decrease values; and the TAB key is used to move the cursor to the next field on a display.

1.5.3 Panel Overview

The Phoenix panel is available in three formats: ID only, AST only, and combination ID/AST testing. The “pour and seal” serpentine design is optimized for safety and leak-resistance. Each well in the disposable contains approximately 50 μ L of inoculum in an environment that prevents significant evaporation during the course of incubation.

A Panel inoculation station holds six tubes of broth (ID, AST) and three panels held at an angle of 24° in order to provide proper gravity-driven inoculum flow through the panel. See Figure 1–2.

A Panel Transportation Caddy is a molded plastic tray used to transport filled and sealed Phoenix panels from the preparation bench to the Phoenix instrument. The caddy capacity is 20 panels. See Figure 1–3.


Figure 1–2 – Panel Inoculation Station


Figure 1-3 – Panel Transportation Caddy

1.5.4 Testing Overview

After a door-opening/-closing event, the instrument reads panel barcode labels and then performs a scan using the red LEDs to determine if panels are present, and to locate or “map” well positions. Panel readings are made on the hour and at 20 and 40 minutes after the hour. The instrument test sequence begins with the system checking that the door is locked. Tiers perform dark readings, then the UV bulb is turned on and allowed to warm up. UV readings are then taken for one revolution. The UV bulb is turned off. Next the red LEDs are turned on, allowed to warm up, and red LED readings are taken. Then the green LEDs go through the same sequence. Finally, the blue LEDs go through this sequence.

1.5.5 Normalizers

Normalizers serve as reference panels for adjustment of the instrument's optical detection system. They are used for adjustment of the LEDs (Red, Green and Blue) and fluorescence lamps.

For Red, Green and Blue correction, data from the normalizer panels are used to correct for variations in optical channel gain and to compensate for well-to-well parallax. Raw signals are source monitor corrected and are then ratioed to the value of the corresponding well of the normalizer panel. This ratio is multiplied by a correction factor, which is the expected transmittance level of the normalizer, to scale the resultant values. The fluorescence (UV) signal from the Normalizer is used to ensure that the UV intensity reaching the panel is within the acceptable range for proper interrogation of the fluorescent ID substrates on the Phoenix panel. If the UV level is either too high or too low at the Normalizer, a UV lamp adjust will be performed at the next time period when there are no ongoing panels. Extreme deviations from the normal UV signal will cause a system alert, resulting in the aborting of panels in the offending tier(s). This is followed by an attempt to perform a UV lamp adjustment.

Normalizers expire two years after installation. A system alert will occur before expiration to allow you to schedule replacement.

The normalizer panels are always located in station number 0 of each tier. When the door is opened, routine carousel access does not present the normalizer panel. The normalizer panel is constrained in its carrier to prevent inadvertent removal.

Automatic Adjustment of Light Sources

Based on the results of the system's Built-In-Test (BIT), as well as time and power cycling factors, the system may detect that readings from the ultraviolet (UV) and visible light sources (testing lamps) are out of tolerance. There are two main levels to this condition:

- 1 The deviation is great enough that panel results are invalidated (panel testing for the tier is aborted); the tier's stations are blocked; an automatic light source adjustment is performed as soon as all panels in the instrument complete testing
- 2 The deviation is within limits that do not affect panel results; however available stations are blocked; as soon as all panels on the tier complete testing or are removed, an automatic light source adjustment is performed

"E" type error codes and subcodes inform you of both of these conditions, for each of the tiers.

Automatic light source adjustment attempts to bring readings within acceptable ranges. It is performed when the following conditions are met: the instrument is warmed up, idle, and there are no ongoing panels. UV adjustment prohibits access to the instrument for entering panels into unaffected tiers, removing panels, performing maintenance checks, etc.

1.6 Manual Structure

This user's manual contains eight sections, a glossary, appendices, an index, and a Reader Comment Card.

Section ❶ – Introduction – provides an overview of the Phoenix system and its uses in the microbiology laboratory. An overview of this manual's structure and conventions is also included, along with a summary of warnings and cautions that appear throughout the manual.

Section ❷ – Installation – gives specifications for installation of the Phoenix instrument and instructions for instrument setup.

Section ❸ – Controls and Indicators – explains the use and meaning of all controls and indicators of the system.

Section ❹ – Operation – explains how to use the Phoenix system for testing panels and obtaining results.

Section ❺ – Reference – provides detailed reference material on all displays, reports, icons, and other elements of the user interface.

Section ❻ – Maintenance – explains all user system maintenance, including parts replacements.

Section ❼ – Troubleshooting – provides a convenient guide identifying error codes and suggesting corrective actions.

Section ❽ – Performance and Limitations – provides tables of performance characteristics and limitations of the procedure.

The Glossary explains several Phoenix system and computer terms used in this manual, as well as abbreviations.

The Appendices contain supplemental information, such as warranty, replacement parts list, a software update log, a listing of BD international contacts, panel information, etc.

The Index provides a listing of major topics and associated page numbers.

1.7 Use of this Manual

This user's manual is designed as an integral part of system operation for technologists, supervisors, and other personnel who operate and maintain the Phoenix instrument on a regular basis. Every attempt has been made to include all information which would be required during normal use and maintenance of the system. Should a question arise which is not answered in this manual, please contact the following parties (USA):

For assistance with mechanical, electrical, or software performance problems, in North America:

☎ Field Service 1-800-544-7434

For procedural, panel, or software operation questions, in North America:

☎ Technical Services 1-800-638-8663

International contacts are listed in Appendix D, or contact your local BD representative.

Comments or recommendations on this user's manual may be expressed on the postage-paid Reader Comment Card at the end of the manual, or send e-mail to TechWriting@bd.com.

Other documentation required for proper system operation includes:

Panel Package Inserts – These documents contain important information on panel use, performance, quality control, and limitations of procedures. A package insert is included with each box of panels, and is available upon request from BD.

1.8 Conventions

1.8.1 Keys

The six hard keys that have fixed functions are the UP ARROW (Increase) key, the DOWN ARROW (Decrease) key, the SILENCE ALARM key, the TAB key, the LOAD PANELS key, and the UNLOAD PANELS key. These keys are marked with symbols representing their functions, and operate identically regardless of the active display or operation. Eight other keys, the software (soft) keys are teal in color, and have functions that vary depending on the active display. Each display shows icons representing the current soft key assignments at the bottom of the screen. To perform the function represented by the icon, press the corresponding soft key, located just below the icon defining it. You can also press the corresponding function key on the keyboard to perform the function (e.g., **F1** corresponds to Soft Key number 1, the leftmost soft key).

The six fixed function keys are always identified in your Phoenix manual by CAPITAL LETTERS (e.g., SILENCE ALARM key). The software keys are always identified by lowercase letters in quotes, and the words *soft key* (e.g., "panel results" soft key). Field names are shown in **bolded text**.

1.8.2 Symbols Used on the Equipment

The following symbols appear on the side of the Phoenix instrument:


Figure 1-4 – Symbols Used on the Phoenix Instrument

Top figure: Clockwise from top left, symbols for Serial Port, Printer Port, Network Port, Remote Alarm Port, Barcode Scanner Port; Bottom figure: Symbols for (left to right): electrical hazard; "refer to accompanying documentation" (specifically, the user's manual) for instructions; On/Off; Biohazard; Protective Conductor Terminal

1.8.3 Notes, Cautions, and Warnings

Throughout this manual, important information is presented in boxes offset from the regular text, and is labeled as either a NOTE, CAUTION, or WARNING. These messages are formatted as shown below and bear the following significance:

NOTE

Important information about system use worthy of special attention is presented as a NOTE.

CAUTION

Information on an activity which potentially could cause damage to the instrument or system is presented as a CAUTION.

WARNING

INFORMATION ON AN ACTIVITY WHICH POTENTIALLY COULD CAUSE INJURY TO THE USER IS PRESENTED AS A WARNING.

1.9 Summary of Cautions and Warnings

- Protection provided by this equipment may be impaired if the equipment is used in a manner not consistent with the instructions in this manual.
- The intake filter on the top of the Phoenix instrument must remain unobstructed at all times. Restricted air flow may cause excessive temperatures in the instrument, which can affect test results and possibly cause hardware malfunctions.
- Due to the size and weight of the Phoenix instrument, never attempt to lift the instrument without the aid of mechanical lifting devices. Serious injury can result from attempting to lift the instrument.
- All system users should become thoroughly familiar with all controls and indicators before attempting to operate the instrument.
- Observe established precautions against microbiological hazards throughout all procedures. All specimens should be handled according to CDC-NIH recommendations, CLSI guidelines, or local institution guidelines for any potentially infectious human serum, blood, or other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.
- In addition to wearing gloves, the use of disposable lab coats or gowns and protective glasses or goggles is recommended when working around the instrument.

- The instrument door is electromechanically latched and is controlled by the instrument software. Never attempt to defeat the door latching mechanism, or to open the door when the "unlocked" icon is not displayed. Serious injury can be caused by the rotating carousel. If the carousel is not completely stopped when the door is opened, immediately contact BD for service. Never attempt to rotate the carousel manually or serious injury may result.
- When the system notifies you of alerts and errors, you should immediately respond to the condition.
- All maintenance and repair other than the procedures described in Section 6.2 – Routine Maintenance and Section 6.3 – Module Replacement, must be performed by qualified service personnel. Non-compliance with this warning may result in personal injury or instrument malfunction.
- If you return the existing filter, insure that you replace the filter in the same orientation that it was removed (i.e., with arrow stamped on the side of the filter pointing in the same direction).
- All portions of the body that could possibly come into contact with the affected instrument surfaces must be completely covered before beginning the decontamination process.
- If any error sub-codes other than those listed here appear, note the sub-code and contact BD for assistance.
- If the recommended corrective actions do not solve the problem, contact BD.

2 Installation

2.1 General

This section provides specifications for installation and setup of the BD Phoenix system. The following major topics are included:

- ◆ Instrument specifications
- ◆ Software setup
- ◆ Software installation

WARNING

PROTECTION PROVIDED BY THIS EQUIPMENT MAY BE IMPAIRED IF THE EQUIPMENT IS USED IN A MANNER NOT CONSISTENT WITH THE INSTRUCTIONS IN THIS MANUAL.

CAUTION

The intake filter on the top of the Phoenix instrument must remain unobstructed at all times. Restricted air flow may cause excessive temperatures in the instrument, which can affect test results and possibly cause hardware malfunctions.

2.1.1 Installation

The BD Phoenix instrument is installed only by BD representatives.

WARNING

DUE TO THE SIZE AND WEIGHT OF THE PHOENIX INSTRUMENT, NEVER ATTEMPT TO LIFT THE INSTRUMENT WITHOUT THE AID OF MECHANICAL LIFTING DEVICES. SERIOUS INJURY CAN RESULT FROM ATTEMPTING TO LIFT THE INSTRUMENT.

The BD Phoenix instrument should be installed in an area that is free from undue vibration, direct sunlight, high humidity, dust, temperature extremes, and corrosive or explosive vapors or gases. The system will operate within specifications in room temperatures between 18° – 30° C (64.4° – 86° F). Relative humidity should be between 20% and 90% (non-condensing). The left, rear, and right sides of the instrument should be placed at least four inches from any wall. Environments which exceed these limits could adversely effect the performance of the system components.

The carousel should maintain its temperature to within plus or minus 1.5° C of the temperature controller's setting (35° C). This accuracy can be assured only if the room temperature meets the requirements given above.

Use of earthquake anchoring is strongly recommended in areas susceptible to earthquake activity.

2.2 Instrument Specifications

Physical Dimensions	
Height	41.7 in (105.9 cm)
Width	44.3 in (112.5 cm)
Depth	29.5 in (74.9 cm)
Clearance – Right, Left, Rear	4 in (10.2 cm)
Clearance – Front	18 in (45.7 cm)
Weight – Empty	333 lb (151.4 kg)
Weight – Full	351 lb (159.6 kg)
Counter must be capable of supporting 500 lb (227.3 kg)	

Table 2-1 – Physical Dimensions

Environmental Requirements	
Non-Operating Storage	
Temperature	-17.8° C —> 65.0° C (0° F —> 149° F)
Humidity	10% – 90% RH, non-condensing
Operating Conditions	
Temperature	18° C – 30° C (64.4° F – 86° F)
Humidity	20% – 90% RH, non-condensing
Locations	Level Surface, No direct sunlight, No direct heat
Use of earthquake anchoring is strongly recommended for the Phoenix instrument in locations susceptible to earthquake activity.	
Installation Category II and Pollution Degree 2 as per IEC 664.	

Table 2-2 – Environmental Requirements

Electrical Requirements	
Input Voltage	100 – 117 VAC or 220 – 240 VAC ±10%
Input Current	8 Amp
Input Line Frequency	47 – 53 Hz or 57 – 63 Hz
Power	1.03 Kwatts
Heat	3,515 Btu/hr.

Table 2-3 – Electrical Requirements

Optical Specifications	
Peak Wavelengths	
428 nm to 623 nm, visible spectrum	
367 nm, UV Excitation	
410 nm -- 640 nm Bandpass, UV Emission	

Table 2-4 – Optical Specifications

2.3 Software Setup

The system ships with all setup parameters preset to factory default values. However, before using the instrument for panel testing, you should review the setup parameters to see if they are suitable for your laboratory. These parameters are grouped in the following categories:

- ◆ Instrument Configuration (Section 2.3.1)
- ◆ Communications Configuration (Section 2.3.2)
- ◆ Custom Interpretation Rule Set (Section 2.3.3)
- ◆ BDxpert Rule Configuration (Section 2.3.4)
- ◆ Rapid Reporting Configuration (Section 2.3.5)
- ◆ Panel Lot Definition (Section 2.3.6)

Instrument setup parameters are explained in the paragraphs below. The Main Status screen, which is active immediately after instrument startup, is shown in Figure 2-1.


Figure 2-1 – Main Status Screen

To access the Configuration Menu, press the soft key corresponding to the icon shown below (the "configuration" soft key).

When you press the "configuration" soft key, the password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key. All users with valid passwords have access to general instrument configuration parameters. Only administrators have access to password administration functions.

The configuration menu selections are presented as a new series of soft keys (Figure 2-2). Press the soft key corresponding to the type of configuration function you want to perform. The soft keys are as follows:

"instrument configuration"

"maintenance"

"communications configuration"

"custom interpretation rule set configuration"

"BDXpert rule configuration"

"rapid reporting configuration"


Figure 2-2 – Configuration Menu

Maintenance functions, accessed through the second soft key, are described in other sections of this manual where appropriate.

"panel lot definition"

"exit"

Any changes to configuration parameters are in effect from the time of the change forward.

* Note that the "custom interpretation rule set configuration" and "BDXpert rule configuration" icons do not appear if the Phoenix instrument is connected to and communicating with an EpiCenter system.

2.3.1 Instrument Configuration

To access general instrument configuration functions, press the “instrument configuration” soft key.

The Instrument Configuration display then appears (Figure 2–3).

When the desired configuration parameters have been entered or modified, press the “save” soft key to save the information.

2.3.1.1 General Instrument Configuration

The following parameters can be set in Instrument Configuration:

GENERAL Window

The GENERAL window is located at the top left of the display.

Instrument Number

Select the instrument identification number. The default setting is 1. To increase or decrease the instrument number, use the UP ARROW or DOWN ARROW key. You can choose a number from 1 to 99. If there is only one instrument at your location, you should leave this value set at 1.

Language

Select the language for screen displays and instrument reports. The default setting is English. To scroll through the available selections, use the UP ARROW or DOWN ARROW key. You can choose from the following language selections (shown in the specified language):

English	English	Español	Spanish
Français	French	Português	Portuguese
Deutsch	German	日本語	Japanese
Italiano	Italian		

Rule Set

Select the rule set you want the interpretation engine to use. Only one rule set can be used. You can choose from the following selections:

Custom (defined in Custom Interpretation Rule Set, Section 2.3.3)

EUCAST (European Committee on Antimicrobial Susceptibility Testing)

CLSI (Clinical and Laboratory Standards Institute)

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Rule Version

This read-only field shows the current version of the rule set selected in the previous field.

(Audible) Alarm Volume

Select the volume of the instrument's audible alarm. The default setting is 5, which is at the middle of the volume range. To increase or decrease the volume, use the UP ARROW or DOWN ARROW key (a sample volume tone sounds each time you adjust the setting). You can choose from 0 (audible alarm off) to 10 (loudest). Only the volume of "Alert" and "Activity" alarms (see **Section 3.7 – Audible Tones and Alarms**) is affected by this setting.

Resistance/Critical Checkbox

This checkbox appears *only* when the Phoenix instrument is connected to a BD EpiCenter system but is *not* communicating with it. Check this box to enable the Phoenix instrument to operate in a standalone mode when communications with EpiCenter is interrupted. When this box is checked, the **BDXpert System Active** checkbox appears.

BDXpert System Active Checkbox

Press the **[Spacebar]** or "select" soft key to check this box. Checking the box activates the entire BDXpert System rules, which comprises CLSI, EUCAST, SFM, or Custom Interpretation Rules.

This field does not appear if the EpiCenter system is attached and communicating. However, if communication with the EpiCenter system is lost, the checkbox reappears. This enables you to reactivate BDXpert rules interpretations in the standalone Phoenix system. When communications with the EpiCenter system is restored, BDXpert rules interpretations are once again performed at the EpiCenter system, and this checkbox is removed from the display.

NOTE

Disabling all BDXpert rules also disables detection of Resistance Markers (e.g., ESBL), except those triggered by 1500-series rules.

INSTRUMENT CONFIGURATION	
GENERAL Instrument # 1 Language English Rule Set CLSI Rule Version Alarm Volume 5 BDXpert System Active ☒ Rapid Completion ☒	DATE/TIME Date Format mm/dd/yy Time Format hh:mm Daylight Saving Time ☐ Time Zone Offset -5 : 00 Current Date 11 / 16 / 07 Current Time 14 : 29
OPTIONS Abbreviated Lab Report ☐ QC Lot Support ☐	Daily Instrument Report Printing Phoenix ☐

Figure 2-3 – Instrument Configuration Display

Rapid Completion Checkbox

Press the **[Spacebar]** or “select” soft key to check this box. Checking the box enables the instrument to provide BDXpert AST results (SIR) before determining actual MIC values. The instrument MIC values are provided as soon as they can be accurately determined. Within a test panel, some MIC values may be available earlier than others. The rapid completion feature can be used to predict resistance for uncompleted antibiotics using the ID alone (intrinsic resistance), or ID with completed MICs of related antibiotics, or resistance marker tests (BL, ESBL). The BDXpert system is used to make these predictions. This can be useful in situations where, for example, the results for drugs that have not yet received MICs would be of no clinical value based on the other results that are already available. Antimicrobials with Rapid Complete BDXpert interpretations are indicated by a “C” in the MIC column on Results screens and Lab reports. Uncheck the checkbox to disable rapid completion.

If Rapid Completion is enabled at the EpiCenter system, then it must also be enabled (checked) at the Phoenix instrument.

OPTIONS Window

The OPTIONS window is located at the bottom left of the display.

Abbreviated Lab Report

Press the **[Spacebar]** or “select” soft key to check this box. Checking the box causes the system to print a shortened version of the Lab Report. The shortened version does not contain results for the biochemical (ID) reactions. Uncheck the box to print the standard full-length Lab Report. The default is unchecked (full length). Note that QC Reports always print standard full-length Lab Reports.

QC Lot Support

Press the **Spacebar** or “select” soft key to check this box. Checking the box enables the QC Lot feature, which can be used to facilitate panel lot QC testing and tracking. When enabled, a new soft key (“panel lot definition”) appears on the Configuration Menu. Refer to Section 2.3.6 **Panel Lot Definition** for more information.

DATE/TIME Window

The DATE/TIME window is located at the top right of the display.

Date Format

To adjust the date format, press the DOWN ARROW key to drop down a listing of all the available date formats. The default format is MM/DD/YY. Use the UP ARROW or DOWN ARROW to highlight the desired selection. You can choose from the following:

Slash separators (/)	MM/DD/YY or MM/DD/YYYY	DD/MM/YY or DD/MM/YYYY	YY/MM/DD or YYYY/MM/DD
Hyphen separators (-)	MM-DD-YY or MM-DD-YYYY	DD-MM-YY or DD-MM-YYYY	YY-MM-DD or YYYY-MM-DD
Period separators (.)	MM.DD.YY or MM.DD.YYYY	DD.MM.YY or DD.MM.YYYY	YY.MM.DD or YYYY.MM.DD

Table 2-5 – Date Formats

Press the **Spacebar** or “select” soft key to select the highlighted entry.

Time Format

To adjust the time format, press the DOWN ARROW key to drop down a listing of all the available time formats. The default format is hh:mm. Use the UP ARROW or DOWN ARROW to highlight the desired selection. You can choose from the following: hh:mm OR hh:mm OR hh:mm.

Press the **Spacebar** or “select” soft key to select the highlighted entry.

Daylight Saving Time Checkbox

Press the **Spacebar** or “select” soft key to check this box. Checking the box adjusts the instrument time for daylight saving time. Uncheck the box to return to standard time. This field does not appear if the EpiCenter system is attached and communicating.

Time Zone Offset

Use the UP ARROW or DOWN ARROW to set your time zone difference from GMT (Greenwich Mean Time). Select negative integers if you are west of the Prime Meridian and east of the International Date Line. Select positive integers if you are east of the Prime Meridian and west of the International Date Line. This field does not appear if the EpiCenter system is attached and communicating.

Current Date

The current date is displayed in the format selected in the Date Format field.

To adjust the month, day, or year, press the “tab” key (or **Tab** key) to highlight the appropriate value according to your date format. Press the UP ARROW or DOWN ARROW key to increase or decrease the displayed value. This field does not appear if the EpiCenter system is attached and communicating.

Current Time

When this display is first accessed, the hours value is highlighted. To adjust the hours, use the UP ARROW or DOWN ARROW key to increase or decrease the displayed value.

To adjust the minutes value, press the “tab” key to highlight the minutes field. Use the UP ARROW or DOWN ARROW key to increase or decrease the displayed value. This field does not appear if the EpiCenter system is attached and communicating.

You cannot change the time while panels are in protocol.

Daily Instrument Report Printing Window

The Daily Instrument Report Printing window is located at the bottom right of the display.

Phoenix Checkbox

Press the **Spacebar** or “select” soft key to check this box. Checking the box causes the Auto Time field to appear, and tells the instrument to print a Daily Instrument Report automatically at the time specified in the **Auto Time** field. Uncheck the box to disable automatic printing. (If you disable automatic printing, you can still print the report at any time from the Reports Menu.)

Auto Time

In the first Auto Time field, press UP ARROW or DOWN ARROW to enter the hour at which you want the Daily Instrument Report to print automatically. Time is calculated in 24-hour international (military) format. Press **Tab** to advance to the next time field and press UP ARROW or DOWN ARROW to enter the minutes.

2.3.1.2 More Configuration Options

A display with additional configuration options is accessible by pressing the “more configuration” soft key from the Instrument Configuration display. See Figure 2–4.

INSTRUMENT CONFIGURATION									
Columbia Agar w/5% SB									
LABORATORY INFORMATION									
Name		Kirk Memorial Hospital							
Address		1313 Kirk Boulevard							
City		Kirkville, USA							
Director		Dr. LJ Kirk							
INOCULUM DENSITY					DEFAULT MEDIA TYPE				
		0.25		0.5					
Gram Negative		☐		☒			Yeast		
Gram Positive		☐		☒			COLSB		
							CHOC		
							COLSB		
							SABDX		
							SABEM		
							SABHI		
							TSASB		
 ☐ ☒ ☐									

Figure 2-4 – More Configuration Options Display

Laboratory Information

Information you enter in the Laboratory Information window prints in the header of system reports.

Name

Enter your laboratory name, up to 40 characters.

Address

Enter your address, up to 40 characters.

City

Enter your city (as well as state and zip code, if desired), up to 40 characters.

Director

Enter the name of the laboratory director, up to 40 characters.

Inoculum Density

The Inoculum Density window enables you to set the default McFarland concentration for inoculum. You may select 0.25 (acceptable density of 0.20 – 0.30) or 0.5 (acceptable density of 0.50– 0.60) for Gram Positive and Gram Negative panels (Strep panels use only 0.5). The default density is 0.5 for Gram Negative and Gram Positive panel types. The density for Yeast ID panels is fixed at 2.0 McFarland (acceptable density of 2.00 –2.40). Inoculum density is applicable only to panels that have an ID side.

To select a different density, tab to the desired density and press the "select" soft key. A check appears in the selected density checkbox.

Default Media Type

Yeast – This field enables you to select a default media type that appears during Panel Login when a Yeast ID panel sequence number is scanned or typed in. Press ▼ to drop down a box listing all media types (abbreviations), sorted alphabetically. (When you highlight the media type, the full name appears at the top right of the display.) Use ▼ or ▲ to highlight the desired media. Press the “select” soft key or  to select that media. You do not have to specify a default media type, but you must select a media type when logging in Yeast ID panels.

2.3.1.3 Password Administration

The Password Administration display enables you to set passwords for the three security levels of access. Only users with Administrator level passwords can access the Password Administration display. See Figure 2-5.

The following security levels for passwords exist:

Administrator – Highest access level: access to all user functions and displays

Supervisor – Intermediate access level: cannot access password administration display; access to all other functions and displays

General – lowest access level: cannot access password administration display; cannot access event log display; access to all other functions and displays

Note that the passwords entered for each level must be unique.

To access Password Administration:

From the Main Status screen, press the “configuration” soft key

The password window appears. To access configuration functions, enter the appropriate password then press the “perform action” soft key.

From the Configuration Menu, press the “instrument configuration” soft key

From the Instrument Configuration display, press the “password administration” soft key

Password Administration soft keys:

“save”

“exit”

PASSWORD ADMINISTRATION	
Administrator Password	CAPTKIRK
New Password	
Supervisor Password	ENTERPRISE
New Password	
General Password	BEAMMEUP
New Password	
 	

Figure 2-5 – Password Administration Display

To assign a new password:

- 1 The Administrator Password shows the current administrator level password (ADMIN by default). To change this password, type the new password in the New Password field, from 5 to 20 alphanumeric characters. The password is **not** case sensitive. Press the “save” soft key to save the changes. The new password immediately appears in the Administrator Password field.
- 2 The Supervisor Password shows the current supervisor level password (LISVIEW by default). To change this password, type the new password in the New Password field, from 5 to 20 alphanumeric characters. The password is **not** case sensitive. Press the “save” soft key to save the changes. The new password immediately appears in the Supervisor Password field.
- 3 The General Password shows the current general level password (ECOLI by default). To change this password, type the new password in the New Password field, from 5 to 20 alphanumeric characters. The password is **not** case sensitive. Press the “save” soft key to save the changes. The new password immediately appears in the General Password field.

2.3.2 Communications Configuration

The Communications Configuration display allows you to enable/disable communications from the Phoenix instrument to a compatible LIS (Laboratory Information System), and to adjust communications parameters for this type of connection. (It also enables BD representatives to enable/disable/adjust communications with the BD EpiCenter advanced data management system.) Only LIS or BD EpiCenter communications can be enabled. If you have a BD EpiCenter system and also need to communicate with a LIS system, BD EpiCenter can be configured to communicate with the LIS.

When the LIS communications state is changed (from disabled to enabled), the Phoenix instrument must be rebooted for the changes to take effect. In addition, changes to any of the LIS communications setup parameters require an instrument reboot to take effect.

LIS operations are discussed in Section 4.17. This section explains configuring LIS communications only. See Figure 2–6.

To access Communications Configuration:

From the Main Status screen, press the “configuration” soft key

The password window appears. To access configuration functions, enter the appropriate password then press the “perform action” soft key.


From the Configuration Menu, press the “communications configuration” soft key

Communications Configuration soft keys:

“save”

“select”

“organism configuration”

“antimicrobial configuration”

“exit”

The following parameters can be set in Communications Configuration:

Communications Configuration window

LIS Enabled checkbox

Press the “select” soft key or **[Spacebar]** or **[Enter]** to check the box and enable LIS communications. If LIS communications is enabled and you wish to disable it, press the “select” soft key or **[Spacebar]** or **[Enter]** to uncheck the box.

Options window

Send Interpretation Results checkbox

This checkbox causes the final SIR values for antimicrobials to be included in the Results record uploaded to the LIS. Its default value is checked (enabled). To disable it, press the “select” soft key or **[Spacebar]** or **[Enter]**.

Unsolicited Queries checkbox

This checkbox enables the Phoenix instrument to request panel information from the LIS if the panel is placed into the instrument *lacking* an organism ID for AST panels or Combination panels with only the AST portion of the panel enabled.

This field's default value is unchecked (disabled), which means the Phoenix instrument will NOT request missing information from the LIS. To enable it, press the "select" soft key or **Spacebar** or **Enter**.

Send When Placed in Instrument checkbox

This checkbox causes the Phoenix instrument to send a Results upload to the LIS when the panel is placed into the instrument.

This field's default value is checked (enabled). To disable it, press the "select" soft key or **Spacebar** or **Enter**.

ASTM Byte Mode Comments checkbox

This checkbox only appears when Japanese is selected for the language.

This checkbox causes the ASTM standard to be followed (when enabled/checked) or not followed (when disabled/unchecked). When the ASTM standard is not followed, escape sequence characters that are usually rejected are instead accepted.

This field's default value is checked (enabled). To disable it, press the "select" soft key or **Spacebar** or **Enter**.

Results Upload Options window

Results upload records consist of the following information:

- ◆ Header record (Delimiter fields, sender name, version number, message date/time)
- ◆ Order record (Accession number, Isolate number, Organism, Test ID, Sequence number, Priority, Report type)
- ◆ Comment record (contains Special Messages and/or BDXpert Rules)
- ◆ Results record
 - ◆ Panel sequence number
 - ◆ Instrument number
 - ◆ Instrument type
 - ◆ Instrument location (Station)
 - ◆ Time to identification (if applicable)
 - ◆ Test start time
 - ◆ Test end time
 - ◆ Test status

Pending – Panel has been logged in but has not been detected by instrument inventory.

Ongoing – Panel has been logged in and has been detected by instrument inventory, but no results are available.

Partial complete – One or more results (MIC or ID) have been generated but not all tests on the panel. BDXpert rules may be pending and the panel could have a Needs Attention condition.


COMMUNICATIONS CONFIGURATION	
COMMUNICATIONS CONFIGURATION LIS Enabled ☒	
OPTIONS ☒ Send Interpretation Results ☒ Send When Placed in Instrument ☐ Unsolicited Queries	
RESULTS UPLOAD OPTIONS ☐ Solicited ☐ Send on Finalization ☒ Send on Completion ☐ Send as Available ☐ Send at Fixed Time	COMMUNICATIONS SETUP Baud 9600 Data Bits 8 Parity None Stop Bits 1 Packed Frames Yes
    	

Figure 2-6 – Communications Configuration Display

Complete – The instrument has completed its readings and all panel test results have been determined, or cannot be interpreted. BDxpert rules may be pending and the panel could have a Needs Attention condition.

Rapid complete – The panel has not yet completed testing but the system has determined results for all drugs

Complete with unignored Needs Attention condition – Panel has completed testing and has active needs attention conditions associated. The panel cannot be finalized until the Needs Attention conditions are resolved or ignored.

Complete with all ignored Needs Attention conditions – Panel has completed testing and has Needs Attention conditions that were ignored by the user. By ignoring the Needs Attention condition, the user can finalize the panel.

Complete QC Pass – Panel was logged in as a QC panel, has completed testing, and has passed the QC criteria.

Complete QC Fail – Panel was logged in as a QC panel, has completed testing, and has failed the QC criteria.

- ◆ Result type
- ◆ Antimicrobial code
- ◆ MIC value
- ◆ S / I / R / No Interp / Error value

- ◆ Resistance marker
- ◆ ID or Final ID
- ◆ Results status (finalized or unfinalized)

Only **ONE** of the Results Upload Options can be checked (enabled).

Solicited checkbox

This checkbox causes the Phoenix instrument to upload Results records **ONLY** when the LIS requests the information. Its default value is unchecked (disabled). To enable it, press the "select" soft key or **[Spacebar]** or **[Enter]**.

QC panels and orphan panels are uploaded only when solicited by the LIS.

Send on Finalization checkbox

This checkbox causes the Phoenix instrument to upload Results records **ONLY** when the panel is finalized. Its default value is unchecked (disabled). To enable it, press the "select" soft key or **[Spacebar]** or **[Enter]**.

Send on Completion checkbox

This checkbox causes the Phoenix instrument to upload Results records **ONLY** when the panel status becomes complete **OR** when a change is made to a complete panel. Its default value is checked (enabled). To disable it, tab to another Results Upload option and enable that option.

Send as Available checkbox

This checkbox causes the Phoenix instrument to upload Results records at the following times: when an ID is determined; when an AST result is determined; when there is a change to a panel record that already has at least partial results. Its default value is unchecked (disabled). To enable it, press the "select" soft key or **[Spacebar]** or **[Enter]**.

Send at Fixed Time checkbox

This checkbox causes the Phoenix instrument to upload Results records for a panel with partial results at a fixed time that you specify. Its default value is unchecked (disabled). To enable it, press the "select" soft key or **[Spacebar]** or **[Enter]**.

When you check the box, two additional fields appear that allow you to specify the time for sending Results records. In the first time field, press ▲ or ▼ to enter the hour at which you want the results to be sent. Time is calculated in 24-hour international (military) format. The hour value defaults to 16 (4 p.m.). Press **[Tab]** to advance to the next time field and press ▲ or ▼ to enter the minutes.

NOTE

Send as Available and Send at Fixed Time results should be considered PRELIMINARY until the panel is complete.

Communications Setup window

Baud

Press the ▲ or ▼ to select the baud rate for serial communications with the LIS system. Available choices are: 2400, 4800, 9600 (default), 14400, 19200, 38400.

Data Bits

Press the ▲ or ▼ to select the number of data bits for serial communications with the LIS system. Available choices are: 7, 8 (default).

Parity

Press the ▲ or ▼ to select the parity for serial communications with the LIS system. Available choices are: None (default), odd, even.

Stop Bits

Press the ▲ or ▼ to select the stop bits for serial communications with the LIS system. Available choices are: 1 (default), 2.

Packed Frames

Press the ▲ or ▼ to select whether packed frames can be used for serial communications with the LIS system. Available choices are: Yes (default), No. Packed frames (Yes) enable the Phoenix system to send multiple records per frame. No indicates that one record per frame is uploaded to the LIS.

2.3.2.1 Organism Configuration

The Organism Configuration display allows you to edit the LIS codes for the organisms in the Phoenix instrument database. LIS codes for organisms must be unique.

The main Organism Configuration display enables you to select the organism (or Test Strain for QC organisms) to be edited. A second display is where you perform the actual edits. Organisms are presented alphabetically by organism short name (refer to Appendix E for a short name/full name list); QC organisms are presented alphanumerically with the ATCC strain number appearing first.

To access Organism Configuration:

From the Main Status screen, press the "configuration" soft key

The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.


From the Configuration Menu, press the "communications configuration" soft key

From the Communications Configuration display, press the "organism configuration" soft key

Organism Configuration soft keys:

To edit an organism LIS code:

In the Organism Configuration display (Figure 2-7a), press ▲ or ▼ to highlight the name in the list and press the "select" soft key, Spacebar or Enter. You can move from the organism box to the QC organism box by pressing the Tab key.

In the Organism Edit display (Figure 2-7b), type the desired LIS code, up to 20 alphanumeric characters. Press the "okay" soft key to save the changes.

ORGANISM CONFIGURATION		
Organism ID	BD Code	LIS Code
Achr. piechaudii	ALCPiE	ALCPiE
Achr. species	ACHRSPE	ACHRSPE
Achr. xylos. ssp denitrificans	ALCDEN	ALCDEN
Achr. xylos. ssp xylosoxidans	ALCXYL	ALCXYL
QC Organism ID	BD Code	LIS Code
10094 S. spiritovorum	BD10094	BD10094
10741 E. faecalis	ATCC10741	ATCC10741
11046 E. coli	BD11046	BD11046
11052 P. aeruginosa	BD11052	BD11052

Figure 2-7a – Organism Configuration Display

ORGANISM CONFIGURATION		
Organism ID	BD Code	LIS Code
> Achr. piechaudii	ALCPIE	ALCPIE
Achr. species	ACHRSPE	ACHRSPE
Achr. xylos. ssp denitrificans	ALCDEN	ALCDEN
Achr. xylos. ssp xylooxidans	ALCXYL	ALCXYL

LIS Code

ALCPIE

OK 

Figure 2-7b – Organism Edit Display

2.3.2.2 Antimicrobial Configuration

The Antimicrobial Configuration display allows you to edit the LIS codes for the antimicrobials in the Phoenix instrument database. LIS codes for antimicrobials must be unique.

The main Antimicrobial Configuration display enables you to select the antimicrobial to be edited. A second display is where you perform the actual edits. Antimicrobials are presented alphabetically by antimicrobial name.

To access Antimicrobial Configuration:

From the Main Status screen, press the "configuration" soft key

The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

From the Configuration Menu, press the "communications configuration" soft key

From the Communications Configuration display, press the "anti-microbial configuration" soft key

Antimicrobial Configuration soft keys:

"save"

"select"

"exit"

To edit an antimicrobial LIS code:

In the Antimicrobial Configuration display (Figure 2-8a), press ▲ or ▼ to highlight the name in the list and press the "select" soft key, **Spacebar** or **Enter**.

In the Antimicrobial Edit display (Figure 2-8b), type the desired LIS code, up to 20 alphanumeric characters. Press the "okay" soft key to save the changes.

ANTIMICROBIAL CONFIGURATION		
Antimicrobial	BD Code	LIS Code
Amikacin	AM	
Amoxicillin	AMX	AMX
Amoxicillin-Clavulanate	AMC	AMC
Amoxicillin-Clavulanate (f)	AXC	AXC
Ampicillin	AM	AM
Ampicillin-Sulbactam	SAM	SAM
Ampicillin-Sulbactam (f)	SXA	SXA
Arbekacin	ARB	ARB
Azithromycin	AZM	AZM
Azlocillin	AZ	AZ
Aztreonam	ATM	ATM

Figure 2-8a – Antimicrobial Configuration Display

ANTIMICROBIAL CONFIGURATION		
Antimicrobial	BD Code	LIS Code
Amikacin	AM	AM
Amoxicillin	AMX	AMX
Amoxicillin-Clavulanate	AMC	AMC
Amoxicillin-Clavulanate (f)	AXC	AXC
Ampicillin	AM	AM
Ampicillin-Sulbactam	SAM	SAM
Ampicillin-Sulbactam (f)	SXA	SXA
Arbekacin	ARB	ARB
Azithromycin	AZM	AZM
Azlocillin	AZ	AZ
Aztreonam	ATM	ATM

LIS Code

OK

Figure 2-8b – Antimicrobial Edit Display

2.3.3 Custom Interpretation Rule Set

Custom Interpretation Rule Set Configuration enables you to tailor interpretation rules to the specific needs or requirements of your laboratory. The current rule set (whether default or already customized) is used as the basis for custom rules. You can print the current custom rule set (if it is selected in Instrument Configuration) by selecting "Interpretation Rule Set Report" from the Reports menu (see Section 5.8.8). Note that this configuration function is not available if the Phoenix instrument is connected to and communicating with a BD EpiCenter system.

Periodically, a software update or Phoenix Update Disk operation may install new breakpoints. When custom breakpoints are changed, an icon appears on the Custom Interpretation Rule Set display that enables you to print a Custom Breakpoint Difference Report. This report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed

Body of Report: Rule Set, Previous Version, Current Version, and columns for Antimicrobial (sort order), Organism Group, Organism Name, Rule Set, Previous and Current Standard and Custom: S(usceptible) value, R(esistant) value, and DTG. Breakpoints are listed in rows of the report.

See Figure 2–10.

Rules (Standards) Updates

The chart below depicts the logic of updates to Rules Sets. The top of the chart shows the effect of rules changes when you have customized a rule set. The bottom of the chart shows the rest of the cases for rules updates.

In the chart, the letters A, B, and C represent the entirety of a rule set (standard). For example, CLSI Rule Version M100_S16 might be indicated by value "A." An update, delivered via the Phoenix Update Disk (PUD), might install CLSI Rule Version M100_S17, which might be indicated by the value "B." Your customization of "A" might be represented by "C." Therefore, if you follow the chart below, the effect of the update is that value "C" (your customized rules based on Version M100_S16) would be retained.

	Old Standard	New Standard	Old Custom	Result after PUD New Custom
If you have custom rules	A	A	B	B
	A	B	C	C
	-	-	A	A
	-	A	A	A
	-	A	B	A
	A	-	B	B
No custom rules	A	A	A	A
	A	B	A	B
	-	A	-	A
	A	A	-	A
	A	B	-	B
	A	-	A	-

Table 2-6 – Rules Updates

To print a Custom Breakpoint Difference Report:

From the Custom Interpretation Rule Set configuration screen, press the “print” soft key

To access Custom Interpretation Rule Set configuration:

From the Main Status screen, press the “configuration” soft key

The password window appears. To access configuration functions, enter the appropriate password then press the “perform action” soft key.

From the Configuration Menu, press the “custom interpretation rule set configuration” soft key

The Custom Interpretation Rule Set display then appears (Figure 2–9a).

When the desired configuration parameters have been entered or modified, press the “save” soft key to save the information.

To use the custom rule set for interpretation, go to the Instrument Configuration display (Section 2.3.1), and advance to the **Rule Set** field. Highlight CUSTOM in the drop down box by pressing the DOWN ARROW, then press  or the “select” soft key to select it.

Custom Interpretation Rule Set soft keys:

“save”		“select”	
“mark for deletion”		“delete marked items”	
“print report”		“exit”	

Custom Interpretation Rule Set fields:

Based On

Select the standard rule set upon which the custom rules are to be based. The standard rule sets are: CLSI, EUCAST, and SFM. Press DOWN ARROW to drop-down a box in which you can highlight and highlight the desired rule set. Press Spacebar or the “select” soft key to select the highlighted set.

NOTE

If you select a rule set other than the current rule set, the currently defined *custom interpretation* rules are deleted.

Rules Window

The window in the middle of the screen shows the current rules in the rule set. This list is read-only and can be navigated with UP ARROW and DOWN ARROW. The following values are shown from left to right (items are explained below under **To modify a rule**):

- Marked for Deletion Indicator (“x” indicates you have marked the rule for deletion)
- Antimicrobial
- DTG (Drug Test Group – see Step 4 below)
- Org (anism) Group
- Organism

S (usceptible Value): $\leq$

R (esistant Value): $\geq$

To add a rule:

- 1 Press the "tab" key to highlight the Rules Window.
- 2 Press the End key (keyboard) to move to the very bottom of the window. A blank line is highlighted.
- 3 Press the "select" soft key. The screen shown in Figure 2–9b appears.
- 4 Select the **Antimicrobial**, DTG, Org. Group, Organism, Susceptible, and Resistant values by pressing the DOWN ARROW to drop down a box listing all available values for the field. Use the UP ARROW or DOWN ARROW to highlight the desired value in the list, then press the "select" soft key or  key.

When all the desired values are selected, press the "ok" soft key to accept the new rule and return to the Rules Window.

- 5 Press the "save" soft key to save the new rule.

To modify a rule:

- 1 Use the UP ARROW or DOWN ARROW to highlight the desired rule in the Rules Window. The antimicrobial name appears in the display title area, on the left.
- 2 Press the "select" soft key. The screen shown in Figure 2–9b appears.
- 3 To change the **Antimicrobial**, press the DOWN ARROW to drop down a box listing all available antimicrobials. Use the UP ARROW or DOWN ARROW to highlight the desired drug in the list, then press the "select" soft key or  key.

CUSTOM INTERPRETATION RULE SET

Based On

CLSI

Antimicrobial	DTG	Org Group	Organism	S	R
AM	A	ENTC_IC	{all}	8.000	16.000
AM	A	ENTERIC_IC	{all}	8.000	32.000
AM	A	HAEINF_IC	{all}	1.000	4.000
AM	A	HAEOTH_IC	{all}	1.000	4.000
AM	A	LISMON_IC	{all}	2.000	-
AM	A	LISTOTH_IC	{all}	2.000	-

Figure 2-9a – Custom Interpretation Rule Set – Screen A

CUSTOM INTERPRETATION RULE SET

Ampicillin

Based On

CLSI

Antimicrobial	DTG	Org Group	Organism	S	R
>AM	A	ENTC_IC	{all}	8.000	16.000
AM	A	ENTERIC_IC	{all}	8.000	32.000
AM	A	HAEINF_IC	{all}	1.000	4.000
AM	A	HAEOTH_IC	{all}	1.000	4.000
AM	A	LISMON_IC	{all}	2.000	-
AM	A	LISTOTH_IC	{all}	2.000	-

Antimicrobial

Ampicillin

Org Group

ENTC_IC

Organism

{all}

Susceptible

8.000

 <= Resistant

16.000

 >=

OK

Figure 2-9b – Custom Interpretation Rule Set – Screen B

- 4 To change the DTG, press the DOWN ARROW to drop down a box listing all available DTGs. Use the UP ARROW or DOWN ARROW to highlight the desired DTG in the list, then press the "select" soft key or  key.

DTGs are derived from the AST Standards and are used to categorize antibiotics into distinct groups. For each drug, the groups are specific to the organism group and the recommended utilization for that drug.

DTGs have no significance in Phoenix alone. They are used in conjunction with the BDxpert rules and are necessary for interface with the BD EpiCenter system.

Generally, these are divided into seven groups shown below. Only the A, B, C and U codes may be reported on the Lab report(s). B and C groups are only reported when they are promoted to A by the BDxpert system.

- A** Always tested and always reported.
- B** Usually tested, but not always reported.
- C** Sometimes tested and not always reported.
- U** Urinary tract specific drug. Tested and reported for source urinary tract.
- O** Other drugs that may be tested. Will not be reported unless changed by user.
- I** Investigational drugs. Drugs not approved for clinical use are never reported.
- N** Not grouped by standard. Will not be reported unless changed by user.

These drug testing codes are used for two purposes in the Phoenix system. First, you can alter the drug testing codes to more closely match the antibiotic formulary and drug utilization guidelines within your institution. The initial codes will be determined by the AST interpretive standard that you have selected. Allowable changes are shown below.

Starting Code Allowable (Recommended) Changes

A	B, C, U, O
B	A, C, U, O
C	A, B, U, O
U	A, B, C, O
O	A, B, C, U
I	None
N	A, B, C, U

The second application of the drug testing codes is for the promotion or suppression of drug results for the chartable report. This application is driven by the rules in the BDxpert system. If the rules alter the drug testing codes, these will be reflected on the Lab Report generated by the Phoenix system. The promotion and suppression actions are shown below. At the EpiCenter level, the drug testing codes will be used to determine which drug results actually appear on the chartable report.

Starting Code	Promotion	Suppression
A	Not applicable	C
B	A	C
C	A	Not applicable
U	Not applicable	C
O	Not allowed	Not allowed
I	Not allowed	Not allowed

Custom Breakpoint Difference Report									
03/05/05 09:11									
Rule Set		CLSI							
Previous Version		M100_S12							
Current Version		M100_S13							
Antimicrobial	Organism Group	Organism Name	Rule Set	Previous			Current		
				S	R	DTG	S	R	DTG
AM	ENTERIC_IC	(all)	CLSI	8.000	32.000	A	2.000	4.000	B
			Custom	4.000	32.000	A	4.000	32.000	A
AM	STRPNE_IC	(all)	CLSI	.	.	.	2.000	-	A
			Custom	2.000	-	A	2.000	-	A
AN	ACIN_IC	(all)	CLSI	16.000	64.000	B	2.000	4.000	A
			Custom	16.000	64.000	B	2.000	4.000	A
.	.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.	.
End of Report									

Figure 2-10 – Custom Breakpoint Difference Report

- To change the **Org(anism) Group**, press the DOWN ARROW to drop down a box listing all available groups. "(All)" indicates that all the organism groups are considered in the rule. Use the UP ARROW or DOWN ARROW to highlight the desired group (or "All organisms") in the list, then press the "select" soft key or **Enter** key. (A list of organism groups is provided in Section 5.8.8.)
- To change or select a specific **Organism**, press the DOWN ARROW to drop down a box listing all available organisms. "(All)" indicates that all the organisms in the specified group are considered in the rule. Use the UP ARROW or DOWN ARROW to highlight the desired organism in the list, then press the "select" soft key or **Enter** key. (Note that if an organism group is selected, the organism list displays only organisms within the selected group.)
- To change the **S** (susceptible Value) (i.e., the concentration at or below which the MIC value for the organism or group will return an interpretation of Susceptible), press the DOWN ARROW to drop down a box listing all available concentrations. Use the UP ARROW or DOWN ARROW to highlight the desired value in the list, then press the "select" soft key or **Enter** key.
- To change the **R** (resistant Value) (i.e., the concentration at or above which the MIC value for the organism or group will return an interpretation of Resistant), press the DOWN ARROW to drop down a box listing all available concentrations. Use the UP ARROW or DOWN ARROW to highlight the desired value in the list, then press the "select" soft key or **Enter** key.
- When the desired values for the current rule have been modified, press the "confirm" soft key to return to the Rules Window to select another rule to modify.

OK

- 10 When the desired rules have been modified, press the "save" soft key to save your modifications to the rule set.

To delete a rule:

- 1 Use the UP ARROW or DOWN ARROW to highlight the desired rule in the Rules Window. The antimicrobial name appears in the display title area, on the left.
- 2 Press the "mark for deletion" soft key. This places an "X" to the right of the antimicrobial abbreviation.

X _____

X _____

- 3 Mark any additional rules you want to delete as described in Step 2.
- 4 To delete all the marked rules, press the "delete marked items" soft key.

- 5 Save the deletions by pressing the "save" soft key.

2.3.4 BDXpert Rule Configuration

BDXpert Rule Configuration allows you to enable or disable individual BDXpert rules, or set them to trigger automatically or manually, according to the specific needs or requirements of your laboratory. You can print the existing BDXpert rule set by selecting "BDXpert Rule Set Database Report" from the Reports menu (see Section 5.8.9). Note that this configuration function is not available if the Phoenix instrument is connected to and communicating with an EpiCenter system.

To access BDXpert Rule configuration:

From the Main Status screen, press the "configuration" soft key

The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.


From the Configuration Menu, press the "BDXpert rule set configuration" soft key

The BDXpert Rule Set Configuration display then appears (Figure 2–10).

When the desired configuration parameters have been entered or modified, press the "save" soft key to save the information.

BDXpert Rule Set soft keys:

BDXpert RULE CONFIGURATION		
Rule Set BD STD		
Only results for ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported for fecal isolates of <i>Salmonella</i> spp. or <i>Shigella</i> spp. No other drugs should be reported.		
Rule #	En/Disabled	Auto/Man
101	Enabled	Automatic
103	Enabled	Automatic
104	Enabled	Automatic
105	Enabled	Manual
106	Enabled	Automatic
110	Enabled	Automatic
111	Enabled	Manual


Figure 2-11 – BDXpert Rule Set Configuration

Note that the "auto/manual" soft key does not appear if the rule is disabled. Based On field appears only if Rule Set is "Custom."

BDXpert Rule Set fields:

Rule Set

The rule set currently being displayed is shown in this field. Values are CLSI, SFM, EUCAST, or Custom. Only rules in the selected set appear in the Rule Selection Window.

Based On

If the rule set is Custom, the rule set on which custom values are based is shown in this field.

Rule Window

The rule window, located at the top of the display, shows the text of the rule selected in the rule selection window at the bottom of the display. GREEN text indicates a rule is enabled and automatic. RED text indicates a rule is disabled (auto/manual status is not applicable). BLUE text indicates a rule is enabled and manual.

Rule Selection Window

The rule selection window, located at the bottom of the display, presents a summary of the status of each of the BDxpert rules and enables you to select rules to modify.

To select a rule for modification, use the UP ARROW or DOWN ARROW to highlight the desired rule in the list.

To enable or disable the rule, press the "enable/disable" soft key to toggle the status. When disabled, a rule is not available for use in the system.

To set the rule to automatic or manual trigger, press the "auto/manual" soft key to toggle the status. When set to "auto," a rule will be triggered automatically by the system. When set to "manual," the rule must be accepted by the user before it is applied to results.

While 1500-series rules can be selected for viewing, they are fixed as enabled and automatic. These settings cannot be changed.

2.3.5 Rapid Reporting Configuration

Rapid Reporting Configuration allows you to specify the criteria for notification of critical panels and resistance markers. For Critical Panels, you can be notified by an automatic Lab Report printing, by an audible alarm, and/or by an upload to the LIS when certain types of results are obtained in non-QC panels. The types of results, which you select, are ID only, partial results, and complete results. For Resistance Markers, you designate the type of notification you want to receive when resistance markers (that are enabled) are triggered. For both, you may choose to disable both printing and audible alarms. In this case, the only notification is that the panel is shown in red text on the Panel Inventory display.

Rapid reporting alarms exist in two states: unacknowledged and acknowledged. Unacknowledged alarm panels continue to display in red on the Panel Inventory display. Acknowledgment is achieved when an automatic Lab Report for the panel prints or when the panel is viewed in Panel Results. Audible alarms continue to sound for unacknowledged alarms until the SILENCE ALARM key is pressed.

To access Rapid Reporting configuration:

From the Main Status screen, press the "configuration" soft key


The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

From the Configuration Menu, press the "rapid reporting configuration" soft key

The Rapid Reporting Configuration display then appears (Figure 2–12).

When the desired configuration parameters have been entered or modified, press the "save" soft key to save the information.

Rapid Reporting soft keys:**Rapid Reporting fields:****Critical Panel Notification Window**

The Critical Panel Notification window lets you set parameters for notification on critical panels. Critical panels are any panels for which you want special notification when results become available. Notification can be audible and/or by immediate printing of a Lab Report (or neither). You can mark panels as critical when you log them in via Panel Login, or subsequently on the Panel Results display.

Auto Print

Press the "select" soft key to check (enable) this field. When enabled, the system prints a Lab Report whenever the selected results parameters occur (ID Only, Partial, or Complete – see below). This field is disabled by default (auto print off). Both Auto Print AND Audible Alarm may be enabled simultaneously.

Upload to LIS

This checkbox appears only if LIS Communications is enabled.

Press the "select" soft key to check (enable) this field. The default value is unchecked (disabled). When enabled, the system uploads results records whenever the selected results parameters occur (ID Only, Partial, or Complete – see below).

This field overrides the results upload values set in Communications Configuration (Section 2.3.2).

Audible Alarm

Press the "select" soft key to check (enable) this field. When enabled, the system sounds an audible alarm whenever the selected results parameters occur (ID Only, Partial, or Complete – see below). This field is enabled by default (audible alarm on). Both Auto Print AND Audible Alarm may be enabled simultaneously.

Tone Select

This field enables you to select the tones sounded by the critical panel audible alarm. It appears only if Audible Alarm is checked (enabled). This field is set to Tone 1 by default. Press the ▲ or ▼ key to change from Tone 1 to Tone 2 to Tone 3. The system sounds a sample tone each time the field is selected or its contents are changed.

RAPID REPORTING CONFIGURATION

CRITICAL PANEL NOTIFICATION

☐ Auto Print

☐ Upload to LIS

☒ Audible Alarm

☐ Tone Select

☐ ID Only

☒ Partial Results

☐ Complete

RESISTANCE MARKER NOTIFICATION

☐ Auto Print

☐ Upload to LIS

☐ All Panels

☒ Critical Panels Only

Resistance Marker	En/Disable
BLACT Beta-lactamase	Enabled
ESBL Extended Spectrum Beta-lactamase	Enabled
HLGR High Level Gentamicin Resistant	Enabled
HLKR High Level Kanamycin Resistant	Enabled
HLMUP High Level Mupiricin Resistant Staphylococcus	Enabled
HLPRSP High Level Penicillin Resistant S. pneumoniae	Enabled
HLSR High Level Streptomycin Resistant	Enabled

Figure 2-12 – Rapid Reporting Configuration Display

NOTE

You must check (enable) one of the following three fields in order to receive any critical panel notifications.

ID Only

Press the “select” soft key to check (enable) this field. When enabled, critical panel notification is performed only when an organism ID is detected on a panel. This field is set to disabled by default (alarm/report on partial results). If this field is enabled, **Partial Results** and **Complete** cannot be enabled.

Partial Results

Press the “select” soft key to check (enable) this field. When enabled, critical panel notification is performed when partial panel results are obtained. This field is set to enabled by default (alarm/report on partial results). If this field is enabled, **ID Only** and **Complete** cannot be enabled.

Complete

Press the “select” soft key to check (enable) this field. When enabled, critical panel notification is performed only when complete panel results are obtained. This field is set to disabled by default (alarm/report on partial results). If this field is enabled, **ID Only** and **Partial Results** cannot be enabled.

Resistance Marker Notification Window

This window does not appear if the EpiCenter system is connected.

Auto Print

Press the "select" soft key to check (enable) this field. When enabled, the system prints a Lab Report whenever a resistance marker is triggered. This field is set to disabled by default.

Upload to LIS

This checkbox appears only if LIS Communications is enabled.

Press the "select" soft key to check (enable) this field. The default value is unchecked (disabled). When enabled, the system uploads results records whenever the selected results parameters occur (ID Only, Partial, or Complete – see below).

This field overrides the results upload values set in Communications Configuration (Section 2.3.2).

Audible Alarm

Press the "select" soft key to check (enable) this field. When enabled, the system sounds an audible alarm whenever a resistance marker is triggered. This field is enabled by default (audible alarm on).

Tone Select

This field enables you to select the tones sounded by the resistance marker audible alarm. This field appears only if Audible Alarm is checked (enabled). This field is set to Tone 1 by default. Press the ▲ or ▼ key to change from Tone 1 to Tone 2 to Tone 3. The system sounds a sample tone each time the field is selected its contents are changed.

NOTE

You must check (enable) one of the following two fields in order to receive any resistance marker notifications.

All Panels

Press the "select" soft key to check (enable) this field. When enabled, resistance marker alarms are generated for all panels in which markers occur. This field is disabled by default. If this field is enabled, **Critical Panels Only** cannot be enabled.

Critical Panels Only

Press the "select" soft key to check (enable) this field. When enabled, resistance marker alarms are generated only for critical panels in which markers occur. This field is enabled by default. If this field is enabled, **All Panels** cannot be enabled.

Resistance Marker

The Resistance Marker window, located at the bottom of the display, allows you to enable and disable notification of individual resistance markers. The window shows the abbreviation and text of the Marker, as well as whether it is currently enabled or disabled. By default, all resistance markers are enabled. Resistance marker alarm notification occurs **only** for enabled markers.

To disable notification of a resistance marker, highlight the desired marker and press the "enable/disable" soft key. The En/Disable field changes from Enabled to Disabled.

2.3.6 Panel Lot Definition

The Panel Lot Definition Display enables you to define panel lots to facilitate QC testing and tracking. Panel lot definition saves the lot number, beginning and ending panel sequence numbers, definition date, and first and last used dates for a lot. Current and historical QC results can be viewed and printed. Note that this configuration function is not available if the Phoenix instrument is connected to and communicating with an EpiCenter system.

The display allows you to:

- Define new panel lots
- Recall existing panel lots
- Review current and historical QC results for panel lots
- Extend expiration dates for panel lots
- Print Panel Lot (current, historical) reports

Users with General and Supervisor level passwords may define new panel lots and recall existing panel lots. Users with Administrator passwords may do these things and may extend the expiration date of existing panel lots.

When a panel lot definition is saved, any Test Strain recommendations that exist are displayed in the Current QC Results window. You can toggle between Current and Historical QC Results windows with the "toggle display" soft key. Current results show a summary of the most recent test for each strain for all instruments, up to 20 entries. Historical results show all tests for a strain for the current instrument only, up to 200 entries. In order to display results from other instruments, Panel Lots must be saved at those instruments and restored to the current instrument (see Sections 6.2.4.5/6, Save/Restore Panel Lots).

Field contents of these windows are explained below.

To access Panel Lot Definition:

From the Main Status screen, press the "configuration" soft key

The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

From the Configuration Menu, press the "panel lot definition" soft key

The Panel Lot Definition display then appears (Figure 2–13).

When the desired configuration parameters have been entered or modified, press the "save" soft key to save the information.

Panel Lot Definition soft keys:

"save"

"clear"

"select"

"panel results"


"print"

"toggle display"

"exit"

PANEL LOT DEFINITION					
PQC01					
Panel Lot #	7500001		Definition Date	07/12/04	
Expiration Date	04/01/06		First Used Date	10/22/03	
Start Sequence #	425200001443		Last Used Date	07/12/04	
End Sequence #	425200001468		Where Defined	1	
Current QC Results					
Test Strain	Test Date	Sequence #	Type	QC Status	
25923 S. aureus	07/12/04	425200001444	ID/AST	REVIEW	
29212 E. faecalis	07/12/04	425200001445	ID/AST	FAIL	
29213 S. aureus	07/12/04	425200001446	ID/AST	ERROR	
51299 E. faecalis	07/12/04	425200001447	ID/AST	PASS	

Figure 2-13 – Panel Lot Definition Display

Panel Lot fields:**Panel Lot #**

This field shows the lot number for the panel. **To enter a new panel lot number, you must scan the lot number barcode (next to the L symbol) on the box.** The panel lot number is encoded in the barcode. The lot number cannot be typed in this field when defining new panel lots.

To recall a panel lot definition, press the DOWN ARROW to drop down the list of all panels, or type the first few numbers of the panel lot number to jump to that portion of the list, or scan the panel lot barcode.

Expiration Date

This field shows the lot's expiration date. This information is encoded in the panel lot barcode and is filled in automatically when the panel lot barcode is scanned. The expiration date cannot be typed in when defining new panel lots. This date can only be changed by scanning a new panel lot barcode provided by the BD Quality department (only if you are logged in as an Administrator).

Start Sequence #

This field shows the starting barcode sequence number for the panel lot. To enter a new Start Sequence number, scan the lowest sequence number in the lot (next to the  symbol on the box). This field can only be completed by scanning the panel sequence barcode; the sequence number cannot be typed in this field.

End Sequence #

This field shows the ending barcode sequence number for the panel lot. To enter a new End Sequence number, scan the highest sequence number in the lot (next to the  symbol on the box). This field can only be completed by scanning the panel sequence barcode; the sequence number cannot be typed in this field.

Definition Date

This field shows the date on which the panel lot was defined. It is completed automatically when the panel lot definition is saved. This is a read-only field that cannot be changed by the user.

First Used Date

This field shows the first date on which a panel from this lot was used. It is completed automatically when the panel is placed in the instrument. This is a read-only field that cannot be changed by the user.

Last Used Date

This field shows the last date on which a panel from this lot was used. It is completed automatically when the panel is placed in the instrument. This is a read-only field that cannot be changed by the user.

Where Defined

This field appears only if the panel lot was defined at a different number instrument than the current one. It shows the instrument number where the panel lot was defined. Panel lot definitions can be saved at one instrument and restored to another instrument through the **Save/Restore Panel Lot Definitions** functions in the **Maintenance Menu** (Section 6.2.4.5/6).

Current QC Results Window

This window appears only after a panel lot is saved. When a panel lot definition is saved, any QC Test Strain recommendations that exist are displayed in the Current QC Results window. Only the most recent QC test is shown for each test strain, up to a maximum of 20 tests. Results for other instruments can be shown in this window by saving those results at that instrument and restoring them at the current one (see **Save/Restore Panel Lot Definitions**, Section 6.2.4.5/6).

The following read-only fields are shown in the Current Results window: Test Strain (sort order by name); Test Date; Sequence Number; (Test) Type (ID, AST, or ID/AST); QC Status (None for uncompleted test; Pass, Fail, Review, or Error for completed tests).

You can review results when the panel is in the instrument database as follows: press the TAB key to move to the Current Results window. Use UP/DOWN ARROW to highlight the desired strain, and press the "panel results" soft key when the desired strain is highlighted.

To print the Current QC Panel Lot Report, press the "print" soft key (with the Current QC Results Window showing). In addition to the fields listed above, the report includes standard report header information (Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument number where the report was printed); Panel Lot number; Panel Type; Expiration Date; Instrument number for each strain tested (blank if instrument is the same); and Tech ID. A sample Current QC Panel Lot Report is shown in Figure 2-14.

You can toggle between Current and Historical QC Results windows with the "toggle display" soft key.

Current QC Panel Lot Report

07/15/04 09:11

Instrument #1

Panel Lot #

7500001

Panel Type

PQC01

Expiration Date

04/01/06

Test Strain	Test Date	Sequence #	Type	QC Status	Instrument #	Tech ID
25923 S. aureus	07/12/04	425200001444	ID/AST	REVIEW	.	ljk
29212 E. faecalis	07/12/04	425200001445	ID/AST	FAIL	.	bee
29213 S. aureus	07/12/04	425200001446	ID/AST	ERROR	.	me
51299 E. faecalis	07/12/04	425200001447	ID/AST	PASS	.	up
.	.	.	.	.	.	.
.	.	.	.	.	.	.
.	.	.	.	.	.	.

End of Report

Figure 2-14 – Sample Current QC Panel Lot Report

Historical QC Results Window

This window appears only after a panel lot is saved. Historical results show all tests for a strain for the current instrument only, up to 200 entries. To select a Test Strain, press the TAB key to advance to the Test Strain field, the press DOWN ARROW to drop down the list of strains. Use UP/DOWN ARROW to highlight the desired strain, then press the "select" soft key. When the desired panel is highlighted, you can jump to the QC Panel Results display, or print a Historical QC Panel Lot Report.

The following read-only fields are shown in the Historical Results window: Test Strain (drop down selection box); Test Date (sort order); Sequence #; (Test) Type (ID, AST, or ID/AST); QC Status (None, Pass, Fail, Review, Error); Tech ID.

To print the Historical QC Panel Lot Report, press the "print" soft key (with the Historical QC Results Window showing). In addition to the fields listed above, the report includes standard report header information (Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument number where the report was printed); Panel Lot number; Panel Type; and Expiration Date. A sample Historical QC Panel Lot Report is shown in Figure 2-15.


You can toggle between Current and Historical QC Results windows with the "toggle display" soft key.

Historical QC Panel Lot Report				
07/15/04 09:11				
Instrument #1				
Panel Lot #	7500001			
Panel Type	PMIC/ID			
Expiration Date	04/01/06			
<u>Test Strain</u> – 25923 S. aureus				
<u>Test Date</u>	<u>Sequence #</u>	<u>Type</u>	<u>QC Status</u>	<u>Tech ID</u>
07/12/04	425350000012	ID/AST	REVIEW	ljk
<u>Test Strain</u> – 29212 E. faecalis				
<u>Test Date</u>	<u>Sequence #</u>	<u>Type</u>	<u>QC Status</u>	<u>Tech ID</u>
07/12/04	425350000014	ID/AST	FAIL	bee
<u>Test Strain</u> – 29213 S. aureus				
<u>Test Date</u>	<u>Sequence #</u>	<u>Type</u>	<u>QC Status</u>	<u>Tech ID</u>
07/12/04	425350000016	ID/AST	ERROR	mc
<u>Test Strain</u> – 51299 E. faecalis				
<u>Test Date</u>	<u>Sequence #</u>	<u>Type</u>	<u>QC Status</u>	<u>Tech ID</u>
07/12/04	425350000018	ID/AST	PASS	up
End of Report				

Figure 2-15 – Sample Historical QC Panel Lot Report

2.4 System Startup

Whenever power is applied to the system, it is initialized, performs self-diagnostics, and reports any problems to the user through alerts. If any files are missing or corrupted which would prevent proper operation of the system, the startup process is aborted. If not, the computer loads the operating system and user interface, and begins the warmup period (indicated by the "Instrument is warming up" icon shown below). Afterwards, the system awaits the initiation of panel testing.

The temperature standard panel should be left in the instrument at least 45 minutes before reading it when the instrument is first powered up.

2.5 Software Installation

Software updates are to be performed by BD representatives.

NOTE

Custom Breakpoints

When the software is updated, or when the Phoenix Update Disk operation is performed, if any of the standard interpretation rule sets have changed, new rule sets are installed into the system database. Since custom rule sets are based on standard rule sets, a custom interpretation rule set can be overwritten with new standard rules. If you use custom breakpoints, be sure to print out the Custom Breakpoint Difference Report (Section 2.3.3) after each system update or Phoenix Update Disk operation.

Controls and Indicators

3.1 General

This section describes the meaning and use of the controls and indicators of the Phoenix instrument. The overall layout of the instrument is shown in Figure 3-1. Individual components are illustrated in figures accompanying the related text.

WARNING

**ALL SYSTEM USERS SHOULD BECOME THOROUGHLY FAMILIAR
WITH ALL CONTROLS AND INDICATORS BEFORE ATTEMPTING
TO OPERATE THE INSTRUMENT.**


Figure 3-1 – Phoenix Instrument Layout

3.2 Power Switch

The system power (On/Off) rocker switch is on the right side of the instrument at the bottom (see Figure 3–2). When in the “O” (Off) position, power is removed from the instrument. When in the “I” (On) position, the switch illuminates green and power is applied to the instrument. Power must be turned On for the incubation and testing modules to work. For normal operation, the power should remain On at all times (except during some maintenance procedures).

Figure 3–2 – On/Off Switch

3.3 Keypad and LCD Display

The keypad and LCD display are located on the front of the instrument, at the top right. The keypad is used to issue commands to the instrument. The LCD display presents setup and status information, as well as the keypad definitions that allow you to perform routine operations. See Figure 3–3.

The controls and indicators of the keypad and LCD display are presented in clockwise order from the top right of the module.

3.3.1 DOWN /UP ARROW Keys

The DOWN and UP ARROW keys are located on the right side of the Keypad/LCD Display. The DOWN ARROW key, represented by a down arrowhead, is used to decrease a displayed value, or to scroll downward in a list. The UP ARROW key, represented by an up arrowhead, is used to increase a displayed value, or to scroll upward in a list.


Figure 3-3 – Keypad and LCD Display

3.3.2 TAB Key

The TAB Key is located below the DOWN and UP ARROW keys. It functions like the **Tab** key on the keyboard. When pressed, it advances the cursor to the next field on the display.

3.3.3 Soft Keys

The eight software (soft) keys are located near the bottom of the Keypad/LCD Display, at the center. These keys are teal colored. None of the keys has a fixed function – the functions of the keys vary depending on the current active display. Each display shows icons representing the current soft key assignments at the bottom of the screen. To perform the function represented by the icon, press the corresponding soft key. The soft keys available on each display are explained in Section 5 – Reference.

3.3.4 UNLOAD PANELS Key

The UNLOAD PANELS Key is used to unload any panels from the instrument that are ready to be removed. When pressed, the instrument rotates the carousel to make the first set of Removable panels (indicated by green LED) accessible when the door is opened. It then unlocks the door so that the panels can be removed. If there are additional panels to be removed, close the door and allow the instrument to reposition the carousel to provide access to additional panels.

3.3.5 LOAD PANELS Key

The LOAD PANELS Key is used to load panels into the instrument for testing. When pressed, the instrument rotates the carousel to make the first set of available stations accessible when the door is opened. It then unlocks the door so that the panels can be inserted. If there are additional panels to be inserted, close the door and allow the instrument to reposition the carousel to provide access to additional stations. Panels can be placed in any station that is unoccupied and has no indicator LEDs lit.

3.3.6 SILENCE ALARM Key

The SILENCE ALARM key is located on the left of the Keypad/LCD Display module. When pressed, it turns off the audible alarm for the current "alert" alarm. It does not silence "door open" type audible alarms, and it does not "clear" the error condition. (See **Section 3.7 – Audible Tones and Alarms**, for an explanation of alarm types and tones.)

3.3.7 System Alert/Door Unlocked Indicator

The Alarm Indicator is located at the top of the Keypad/LCD Display, above the display area. This light illuminates or flashes yellow to notify you of certain instrument conditions.

- When the indicator is FLASHING, an alert condition exists. The indicator continues flashing until the error condition is corrected.
- When the indicator is ON, no alert exists and you **can** open the instrument door.
- When the indicator is OFF, no alert condition exists but you **cannot** open the instrument door.

When an alarm occurs, an icon (resembling the one on the Alarm Indicator) appears in the soft key assignments area of the Main Status Screen. Press the soft key corresponding to the System Alert icon. The first error code is shown on the display. Correct any error conditions as soon as possible by following the directions in **Section 7 – Troubleshooting**. You can scroll through the error log with the UP and DOWN ARROW keys.

If there is an audible alarm sounding, you can silence it by pressing the SILENCE ALARM key.

3.3.8 Display Area

The Display area is located at the center of the Keypad/LCD Display. It is used to present information to you, and to show the Soft Key definitions that allow you to perform routine operations. When the instrument starts up, the Main Status Screen is displayed. Other displays appear as you perform various operations.

The LCD Display is programmed to automatically dim after 10 minutes of inactivity. To return the brightness to normal, press any keyboard or soft key.

More information on displays is presented in **Section 5 – Reference**.

3.4 Barcode Scanners

The built-in barcode scanner is located below the Keypad/LCD Display module. The scanner remains on at all times ready to read a barcode. See Figure 3-1.

To scan a barcode, place the panel in the recessed area below the scanner. If necessary, slowly tilt the panel until the acknowledgment beep sounds (indicating that the barcode was scanned successfully).

A hand-operated scanner is also provided with the system to facilitate certain scanning scenarios, such as scanning numerous panels. To scan a barcode with the handheld scanner, press the trigger button on the bottom side of the scanner. A single beep indicates a successful scan. You must press the button for each panel or other barcode you want to scan.

3.5 Floppy Disk Drive

The floppy disk drive is located on the front of the instrument, below the built-in barcode scanner. Its primary purposes are to enable you to save data to floppy disk, and to perform software updates when they are released. See Figure 3-4.

3.5.1 Floppy Disk Indicator

The Floppy Disk Drive Indicator light is on the left side of the drive below the insertion slot. When off it indicates that no activity is occurring in the drive. When on or flashing, it indicates that the disk drive is accessing a floppy disk. **Do not attempt to eject a floppy disk while this indicator is lit.**

3.5.2 Floppy Disk Eject Button

The floppy disk drive eject button is located to the lower right of the insertion slot. When a floppy disk is inserted fully into the slot, this button extends itself. To remove a disk, fully depress the eject button. **Do not attempt to eject a floppy disk while the floppy disk indicator is lit.**


Figure 3-4 – Floppy Disk Drive Controls and Indicators

3.6 Station Indicators

Each accessible station has a set of LED indicators that inform you of the station or panel's status. The status indicators are located at the upper left of the station. See Figure 3-5. The color (red, green, or amber) and state (on or off) indicate the conditions shown in the table below for a given station.

Color	State	Meaning
Amber	ON	Panel found using "locate panel" soft key
Red	ON	Blocked Station (also Temperature Standard Panel)
Green	ON	Removable Panel
All colors	OFF	Ongoing Panel (if panel is in station) OR Available Station (if station is empty)

Table 3-1 – Station Indicators

Figure 3-5 – Station Status Indicators

3.7 Audible Tones and Alarms

Numerous different sounds are generated by the Phoenix instrument as you perform operations. Each of the sounds is unique. These tones are designed to keep you informed about various operational states of the instrument.

Type	Example	Sound
Informational		
Acknowledge	Scanning a barcode	Single short high beep
Sample audible alert	In Configuration mode, the audible alert volume was adjusted	Three tones progressing from high to low ("Figaro")
Carousel halted	The carousel has stopped, the door can be opened	Three tones progressing from low to high
Door closed	The door has been completely closed and latched	Two short high beeps
Activity completed	An operation has been completed	Three short high beeps
Swap floppy	Current floppy disk is full, insert another with free space	Long medium beep
Alarm		
Activity error (volume configurable)	Incorrect barcode scanned	Short high beep then short low beep – sequence repeated four times
System alert (volume configurable)	Optical failure	Single medium beep –one second on, three seconds off, repeating
Error tone	The attempted action cannot be performed	Short high beep then short low beep
Door open alarm	Door has remained open longer than 5 minutes	Continuous shrill trill
Cannot get into instrument	A higher priority activity prevents user access to the instrument	Single low beep
Critical Panel/ Resistance Marker notification	Critical Panel: partial results, panel completes, or ID is determined; Resistance Marker: marker has been triggered	Tone 1 (Medium frequency): Long beep, short beep, long beep, short beep; Tones 2 (low frequency) and 3 (very high frequency): Long beep, short beep, long beep, long beep
Carousel jammed	Carousel motion is impeded	High pitched trill

Table 3-2 – Audible Tones/Alarms

3.8 Keyboard

The keyboard enables you to enter panel demographic information (e.g., accession numbers), as well as performing all the soft key and keypad functions. In addition to the standard alphanumeric key functions, the following key functions exist:

 – Moves cursor to next field; equivalent to TAB (screened) key

  – Moves cursor to previous field

 – Moves up one “page” in a drop-down list field

 – Moves down one “page” in a drop-down list field

 – Toggles the check on/off in a checkbox; inserts a space character in a text field

 – Selects an item in a drop-down box

 – Moves up in a drop-down list field; displays previous selection in a some fields with predefined values (e.g., instrument number); equivalent to UP ARROW screened key

 – Moves down in a drop-down list field; displays next selection in a some fields with predefined values (e.g., instrument number); equivalent to DOWN ARROW screened key

 –  – Corresponds to the eight teal soft keys from left to right

 – Corresponds to LOAD PANELS key

 – Corresponds to UNLOAD PANELS key

 – Corresponds to SILENCE ALARM key

3.9 Printer

For an explanation of controls and indicators on the printer, refer to the manufacturer’s operating instructions furnished separately.

4 Operation

4.1 General

This section describes the routine operation of the Phoenix 100 instrument. The following major topics are discussed:

- ◆ System Maintenance
- ◆ Storage and Handling
- ◆ Preparing Panels
- ◆ Logging in Panels
- ◆ Inserting Panels in the Instrument
- ◆ Panel Testing
- ◆ Responding to Needs Attention Panels
- ◆ Automatic Association of Panels
- ◆ Obtaining Results
- ◆ Panel Finalization
- ◆ Responding to Alarms and Errors
- ◆ Printing Reports
- ◆ Unloading and Discarding Panels
- ◆ Installing New Panel Formats
- ◆ Phoenix Update Disk (PUD)
- ◆ LIS Operations

4.2 Daily System Maintenance

Each day the following maintenance procedures should be performed. The best time to perform maintenance is first thing in the morning, but it may be done any time you find convenient. The following items should be checked:

- 1 Check the printer's paper supply. If the paper supply is low or exhausted, replace the paper as explained in the manufacturer's operating instructions.
- 2 Record the temperature readout shown on the LCD display and that of the temperature standard panel. Note that the temperature panel can be brought into view by selecting one of the instrument LED check functions on the Maintenance menu (see Section 6.2.2). This temperature should be $35^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$.

4.3 Storage and Handling

Phoenix Panels: Panels are individually pouched and packaged in a box of 25. Panels must be stored unopened at room temperature ($15^{\circ} - 25^{\circ}\text{C}$). Do not refrigerate or freeze. Visually inspect the packaging for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Do not use the panel if there is no desiccant or if the desiccant pouch is torn. If stored as recommended, the panels will retain expected reactivity until the date of expiration.

Phoenix ID Broth: Tubes are packaged as 100 tube packs. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store Phoenix ID Broth tubes at $2^{\circ} - 25^{\circ}\text{C}$. Expiration dating is shown on the tube label.

Phoenix AST Broth and Phoenix AST-S Broth: Tubes are packaged as 100 tube packs. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store Phoenix AST Broth tubes at $2^{\circ} - 25^{\circ}\text{C}$. Expiration dating is shown on the tube label.

Phoenix AST Indicator Solution and Phoenix AST-S Indicator Solution: The indicator solution is individually pouched and packaged as a package of 10 dropper bottles. Visually inspect the bottle for cracks, leaks, etc. Do not use if there appears to be a leak, bottle or cap damage, or any change from a dark blue color. **Store the Phoenix AST Indicator Solution at $2^{\circ} - 8^{\circ}\text{C}$.** Each bottle contains enough solution to test up to 100 panels. Expiration dating is shown on the box, pouch, and bottle label and is for unopened bottles. An opened bottle will be stable for up to 14 days if stored at $2^{\circ} - 8^{\circ}\text{C}$. **Be sure the bottle is held vertically when dispensing the AST Indicator Solution**

4.4 Preparing Panels

The Phoenix system is not for use directly with clinical specimens. Only pure culture isolates of aerobic and/or facultatively anaerobic gram-negative, gram-positive, and yeast organisms are acceptable for testing. The test isolate **must** be a pure culture. It is recommended that cultures be 18 to 24 hours old for gram-negative and gram-positive organisms and 18 to 48 hours old for yeast organisms. For AST testing in the Phoenix system, isolates recovered from non-selective media are recommended. It is recommended that media containing antibiotics not be used for organisms to be tested in the Phoenix system except those specifically claimed in the table below. Selective media may inhibit some strains of bacteria and yeast, therefore caution must be used when selecting isolated colonies from these media. Use isolates from a blood agar plate such as Trypticase Soy Agar with 5% Sheep Blood. The Media selected during login of Yeast ID panels refers to the media on which the organism was grown.

Other recommended media are included in the table below:

Recommended Media	Approved Use			
	ID ¹	AST	Strep	Yeast ID
BBL™ CHROMagar™ Orientation	Yes	Yes ¹	X	X
Bromthymol Blue (BTB) Lactose Agar	Yes ⁴	Yes	X	X
Chocolate Agar	Yes	Yes	Yes ²	Yes
Columbia Agar with 5% Horse Blood	Yes	Yes	Yes ³	X
Columbia Agar with 5% Sheep Blood	Yes	Yes	Yes	Yes
Columbia CNA Agar with 5% Sheep Blood (Gram Positives)	Yes	X	Yes	X
Cystine-Lactose-Electrolyte-Deficient (CLED) Agar	Yes ⁵	Yes	X	X
Dey/Engley (D/E) Neutralizing Agar (Gram Negatives)	Yes	X	X	X
Eosin Methylene Blue (Gram Negatives)	Yes	Yes	X	X
Hektoen Enteric Agar (Gram Negatives)	Yes	X	X	X
MacConkey Agar (Gram Negatives)	Yes	Yes	X	X
Phenylethyl Alcohol Agar (Gram Positives)	Yes	X	Yes	X
Sabouraud Brain Heart Infusion Agar - SAB HI (Yeast)	X	X	X	Yes
Sabouraud Dextrose Agar (Yeast)	X	X	X	Yes
Sabouraud Dextrose Agar-Emmons (Yeast)	X	X	X	Yes
Trypticase Soy Agar with 5% Sheep Blood	Yes	Yes	Yes	Yes
Trypticase Soy Agar with Lecithin and Tween 80	Yes	X	X	X
Trypticase Soy Agar without Blood	Yes	X	X	X
Xylose Lysine Desoxycholate Agar (Gram Negatives)	Yes	X	X	X

¹ The use of CHROMagar Orientation may produce false susceptibility results when testing erythromycin with Gram positive organisms. Antimicrobial susceptibility test results should be confirmed using Trypticase Soy Agar with 5% Sheep Blood.

² This media type should *not* be used for Streptococcal identification with SMIC/ID panels. Chocolate Agar may be used for Streptococcal susceptibility testing only.

³ The use of Columbia Agar with 5% Horse Blood may produce significantly higher MIC for SXT with *Streptococcus* species, which may result in false resistance. Antimicrobial susceptibility test results should be confirmed using Trypticase Soy Agar with 5% Sheep Blood.

⁴ The use of Bromthymol Blue Lactose Agar with Gram Positive organisms should be restricted to Staphylococci for both the 0.5 and 0.25 GP systems.

⁵ The use of Cystine-Lactose-Electrolyte-Deficient Agar with Gram Positive organisms should be restricted to Staphylococci for the 0.25 GP system.

Table 4-1 – Recommended Media

Do not use media containing antibiotics for organisms to be tested in the Phoenix system except those specifically claimed in the table above. It is recommended that the applicator swabs are sterile cotton swabs. Polyester swabs are not recommended. The quality of applicator swabs may vary from vendor to vendor and on occasion, loose fibers may dislodge from the swab affecting McFarland readings.

The usefulness of the Phoenix system or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the *Manual of Clinical Microbiology*¹⁷ for specimen collection, transport, and placement on primary isolation media.

Due to variations in inoculum concentrations prepared with McFarland standards, use of the BBL™ CrystalSpec™ nephelometer, the BD PhoenixSpec™ nephelometer, or the BD Phoenix AP™ instrument is required for adjusting the test inoculum prior to use in the Phoenix system.

Instructions for an optional purity check are provided at the end of this subsection.

WARNINGS

OBSERVE ESTABLISHED PRECAUTIONS AGAINST MICROBIOLOGICAL HAZARDS THROUGHOUT ALL PROCEDURES. ALL SPECIMENS SHOULD BE HANDLED ACCORDING TO CDC-NIH RECOMMENDATIONS, CLSI GUIDELINES, OR LOCAL INSTITUTION GUIDELINES FOR ANY POTENTIALLY INFECTIOUS HUMAN SERUM, BLOOD, OR OTHER BODY FLUIDS. PRIOR TO DISCARDING, STERILIZE SPECIMEN CONTAINERS AND OTHER CONTAMINATED MATERIALS BY AUTOCLAVING.

IN ADDITION TO WEARING GLOVES, THE USE OF DISPOSABLE LAB COATS OR GOWNS AND PROTECTIVE GLASSES OR GOGGLES IS RECOMMENDED WHEN WORKING AROUND THE INSTRUMENT.

Materials Required:

- Phoenix Panels
- Phoenix ID Broth or Phoenix Inoculum Broth
- Phoenix AST and/or Phoenix AST-S Broth
- Phoenix AST Indicator Solution and/or Phoenix AST-S Indicator Solution
- Phoenix Panel closures
- Phoenix Inoculation Station
- Phoenix Transport Caddy
- BBL™ CrystalSpec™ nephelometer, the BD PhoenixSpec nephelometer, or the BD Phoenix AP™ instrument
- 25 µL pipettor and tips

Materials Required but Not Provided:

- Gram Stain Reagents
- Sterile Cotton Swabs, Inoculation Loops or Needles
- Nonselective Culture Plated Media (see ¶11, Section 4.4)
- Incubators
- Biohazard Disposable Container
- Markers, etc.
- Vortex mixer

NOTES

Care should be exercised in handling Phoenix panels. You should handle panels by the **sides only** to avoid marking, smudging, or obscuring the bottom or top of the panel in any way.

Accession bar code labels affixed to a Phoenix panel:

- **Must not be of fluorescent material.**
- Should not cover any Phoenix panel reaction wells.
- Should not cover the Phoenix sequence number (panel) barcode.

The procedure that follows describes all the steps in preparing a combination panel for both identification and susceptibility testing. If you are using a combination panel for only ID or only AST testing, note that certain steps are not applicable in the procedure.

General Panel Preparation

If you are using the BD Phoenix AP™ instrument, refer to the *BD Phoenix AP Instrument User's Manual* for panel preparation.

Phoenix Strep panels, Phoenix Yeast ID panels, Phoenix Plus panels, Phoenix Inoculum Broth, and Phoenix MIC panels have separate instructions that appear after Step 19 of these instructions.

- 1 Confirm the Gram stain reaction of the isolate before proceeding with the inoculum preparation for use in the Phoenix instrument. Once the Gram stain reaction is confirmed, select the appropriate Phoenix panel for inoculation.
- 2 Examine the pouch, and do not use the panel if the pouch is punctured or opened. Remove the panel from the pouch. Discard the desiccant. Do not use the panel if there is no desiccant or if the desiccant pouch is torn.

NOTE

Panels must be inoculated within 2 hours of being removed from the pouch.

- 3 Place the panel on the Inoculation Station with the inoculation ports on top and the pad on the bottom.
- 4 Label a Phoenix ID broth tube with the patient's specimen number. Using aseptic technique, pick colonies of the same morphology with the tip of a sterile cotton swab (do not use a polyester swab) or a wooden applicator stick from one of the recommended media.
- 5 Suspend the colonies in the Phoenix ID broth (4.5 mL).
- 6 Cap the tube and vortex for five seconds.


Figure 4-1 – Panel and Closures

- 7 Allow approximately ten seconds for air bubbles to surface. You can tap the tube gently to aid in eliminating bubbles.
- 8 Insert the tube into the CrystalSpec or BD PhoenixSpec nephelometer. Make sure the tube is inserted as far as it will go. (Refer to the BBL CrystalSpec or BD PhoenixSpec nephelometer product insert for correct usage instructions.)
- 9 If the inoculum density is set to 0.5 McFarland for the panel type being run, then a range of 0.50 – 0.60 is acceptable. If the inoculum density is set to 0.25 McFarland for the panel type being run, then a range of 0.20 – 0.30 is acceptable. **If the density of organisms is low**, you can add colonies from the isolate. Re-vortex the sample and reread to confirm that the correct McFarland has been achieved. **If the density of organisms exceeds 0.6 McFarland**, follow the steps below to dilute the broth. It is very important to accurately indicate the level of the liquid in the tube since this volume is needed to adequately fill the wells in the panel.

- a Using a marker, mark the broth level in the over-inoculated Phoenix ID Broth tube.
- b Using a sterile pipette, aseptically add fresh Phoenix ID Broth to the inoculum. Only Phoenix ID Broth may be used to dilute the inoculum.
- c Vortex the tube and allow to sit for 10 seconds.
- d Place the tube in the nephelometer and remeasure the turbidity of the suspension.
 - If the reading is greater than 0.6, repeat Steps b-d.
 - If the reading is 0.5 - 0.6, go to Step e.
- e Using a sterile pipette, aseptically remove excess broth to the original level indicated by the mark on the tube created in Step a.
Remove excess broth to avoid overfilling the panel. Also, do not remove too much broth, as there may be insufficient broth to adequately fill the panel.
- f Broth may now be used to inoculate the Phoenix AST Broth and/or the Phoenix Panel.

NOTES

- Yeast ID panels must be inoculated using a 2.00 - 2.40 McFarland inoculum density.
- Confirm current instrument settings for inoculum density before inoculating panels.
- See instructions below, **ID Inoculum Density Flexibility**, for information on using alternate densities.
- Only the BD PhoenixSpec Nephelometer and BD Phoenix AP instrument can be used to make inoculum densities of 0.25 McFarland
- Standardized bacterial suspension in ID Broth or Inoculum Broth must be used within 60 minutes of preparation.

- 10 If you are performing identification only, proceed to Step 15 and continue the procedure. If you are inoculating a Phoenix Plus Panel, refer to the section below, **Phoenix Plus Panels**.
- 11 Label a Phoenix AST broth tube (8.0 mL) with the patient's specimen number. Add one free-falling drop of AST Indicator solution to the AST broth tube. Invert to mix. DO NOT VORTEX.

NOTES

- Allow AST Indicator Solution to warm to room temperature before dispensing into AST broth.
- The unused portion of the indicator should be returned to 2° – 8° C as soon as possible. **Do not store at room temperature for more than 2 hours. Opened bottles should be discarded after 14 days from initial opening.**
- If volume other than one drop is added inadvertently, discard the tube and use a fresh tube of AST broth.
- After the addition of the Indicator to AST broth, the mixed solution can be stored in the dark, at room temperature, for as long as 8 hours.
 - Tubes must be used within 2 hours after the addition of AST Indicator Solution if exposed to light.

- 12 If an inoculum density of 0.50 – 0.60 was used, transfer 25 μ L of the bacterial suspension from the ID tube into the AST broth tube. If an inoculum density of 0.20 – 0.30 was used, transfer 50 μ L (use 2 shots if utilizing a 25 μ L pipettor) of the bacterial suspension from the ID tube into the AST broth tube.

NOTE

Panels must be inoculated within 30 minutes of the time that the AST broth inoculum is prepared.

- 13 Cap the AST tube and invert several times to mix.
- 14 Wait a few seconds for air bubbles to surface. You can tap the tube gently to aid in eliminating bubbles.
- 15 Pour the ID tube inoculum into the fill port on ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel. If you are using an AST (only) panel, DO NOT inoculate the ID side of the panel. Retain the ID tube for an optional purity check (see below).
- 16 Pour the AST broth inoculum into the fill port on AST side of the panel (85-well side). Allow the fluid to traverse down the tracks before moving the panel.
- 17 Before placing panel closures check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present remove the droplet with absorbant material. The used absorbant material must be decontaminated before discarding.
- 18 Snap on the panel closures. **Make sure that the closures are fully seated. Use 2 closures regardless of panel type.**
- 19 Visually inspect panels to be sure each of the wells is full. Look at both sides of the panel. Make certain that the wells are not overfilled. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

- Panels must be loaded into the instrument within 30 minutes of inoculation (see Section 4.6).
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
 - Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK

It is highly recommended that the purity of both ID and AST inocula be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24-48 h at 35° C under appropriate conditions.

Phoenix Strep Panels

Phoenix Strep panels are for the identification and antimicrobial susceptibility testing of most *Streptococcus* species. Although *Streptococcus* species may be identified in the Gram positive panels, antimicrobial susceptibility cannot be reported when using these panels. The Phoenix Strep panels, which must be used with Phoenix AST-S Broth and Phoenix AST-S Indicator Solution, provide the conditions required for rapid AST testing of most *Streptococcus* species.

- 1 Follow Steps 1 – 9 (General Panel Preparation) to prepare the suspension of bacteria.
- 2 Add one drop of the Phoenix AST-S Indicator to each AST-S broth tube. Invert to mix. DO NOT VORTEX.

NOTES

- AST-S Broth and AST-S Indicator Solution are for use with the Phoenix Strep panels (SMIC/ID, SMIC) only. These reagents are not interchangeable with the AST Broth and AST Indicator Solution used with Phoenix Gram positive and Gram negative panels.
- Allow AST-S Indicator Solution to warm to room temperature before dispensing into AST-S broth.
- The unused portion of the indicator should be returned to 2° – 8° C as soon as possible. Do not store at room temperature for more than 2 hours. Opened bottles should be discarded after 14 days from initial opening.
- If volume other than one drop is added inadvertently, discard the tube and use a fresh tube of AST-S broth.
- After the addition of the Indicator to AST-S broth, the mixed solution can be stored in the dark, at room temperature, for as long as 8 hours.
 - Tubes must be used within 2 hours after the addition of AST-S Indicator Solution if exposed to light.

- 3 Using a pipettor, transfer 25 µL of the standardized bacterial suspension from the ID tube into the AST-S broth tube.

NOTE

Panels must be inoculated within 30 minutes of the time that the AST-S broth inoculum is prepared.

- 4 Cap the AST-S tube and invert several times to mix.
- 5 Wait a few seconds for air bubbles to surface. You can tap the tube gently to aid in eliminating bubbles.
- 6 Pour the ID tube inoculum into the fill port on ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel. If you are using a Phoenix Strep MIC only panel, DO NOT inoculate the ID side of the panel. Retain the ID tube for an optional purity check (see below).
- 7 Pour the AST-S broth inoculum into the fill port on AST side of the panel (85-well side). Allow the fluid to traverse down the tracks before moving the panel.
- 8 Before placing panel closures, check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present, remove the droplet with absorbant material. The used absorbant material must be discarded with biohazard waste.
- 9 Snap on the panel closures. **Make sure that the closures are fully seated.**
- 10 Visually inspect panels to be sure each of the wells is full. Look at both sides of the panel. Make certain that the wells are not overfilled. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

- Panels must be loaded into the instrument within 30 minutes of inoculation (see Section 4.6).
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
 - Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK


It is highly recommended that the purity of both ID and AST-S inocula be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24-48 h at 35° C under appropriate conditions.

ID Inoculum Density Flexibility

You may run the ID portion of a panel in the opposite mode from what is configured by darkening well A-17 on the back of a panel before placing the panel in the instrument. This allows you to run a panel at an inoculum density of 0.20 – 0.30 even if you are configured for a density of 0.5 for that particular panel type. Likewise, you can run a panel at an inoculum density of 0.50 – 0.60 if you are configured for a density of 0.25.

There is no way to alter the density setting during Panel Login. To use a panel in the opposite density mode, using a black Sharpie™ (permanent marker) blacken the **entire A-17** well as shown in the figure below.

Inoculum densities sent from the Phoenix AP instrument (via the EpiCenter system) cannot be changed by darkening well A-17. Inoculum densities sent from the Phoenix AP instrument are invalid if they are received after the completion of the first test cycle, or if the ID portion of the panel is disabled, or if the panel type does not support the inoculum density. If the instrument receives an invalid inoculum density for an ID panel that differs from the panel's determined inoculum density, the ID side of the panel fails, a Needs Attention code is triggered (Can Not Determine Organism ID), and a Special Message is triggered.

The following table outlines how to run the opposite mode from the density for which the instrument is configured.

Current Instrument Inoculum Density Configuration	Inoculum Concentration Desired for Test Panel	Amount of ID Inoculum to Add to AST Broth**	Well A-17
0.50	0.25	50 µL	Blackened
0.25	0.50	25 µL	Blackened
** If also running AST			

Phoenix Yeast ID Panels

Phoenix Yeast ID panels are for the identification of most clinically relevant yeast and yeast-like species.

- 1 Follow Steps 1 – 8 (**General Panel Preparation**) to prepare the suspension of yeast.
- 2 The inoculum density is set to 2.0 McFarland for Yeast ID panels, with a range of 2.00 - 2.40 as acceptable. If the density of organism is low or the density of organisms exceeds 2.40 McFarland, follow the steps discussed in Step 9 (General Panel Preparation) to obtain the correct McFarland density.
- 3 Pour the ID tube inoculum into the fill port on the ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel.
- 4 Before placing panel closures, check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present remove the droplet with absorbant material. The used absorbant material must be discarded with biohazard waste.
- 5 Snap on the panel closures. **Make sure that the closures are fully seated. Use 2 closures regardless of panel type.**
- 6 Visually inspect panels to be sure each of the wells is full. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

- Yeast ID panels must be inoculated using a 2.00 - 2.40 McFarland inoculum density.
- Standardized bacterial suspension in ID Broth or Inoculum Broth must be used within 60 minutes of preparation.

NOTES

- Panels must be loaded into the instrument within 30 minutes of inoculation (see Section 4.6).
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
 - Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK

It is highly recommended that the purity of ID inoculum be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24-48 h at 35° C under appropriate conditions.

Phoenix Plus Panels

Phoenix Plus panels are designed to perform susceptibility testing on an expanded number of antimicrobial agents. To accomplish these susceptibilities, antimicrobial agents are present on both sides of the Phoenix panel. These panels do not have the ability to perform bacterial identification. Because of the design, the inoculation technique is unique and is outlined below. Two tubes of AST broth will be required.

- 1 Follow Steps 1 – 9 (**General Panel Preparation**) to prepare the suspension of bacteria.
- 2 Add one drop of the Phoenix AST Indicator to each AST broth tube.
- 3 Transfer 25 μ L (50 μ L if low inoculum option is used) of the suspension to two Phoenix AST broth tubes. Cap and gently invert.
- 4 Using sterile technique, remove 3.5 mL of broth from one of the inoculated Phoenix AST broth tubes and discard in an appropriate container.
- 5 Pour the remaining 4.5 mL into the left side of the Phoenix Plus panel. Pour the other Phoenix AST broth tube into the right side of the Phoenix Plus panel.
- 6 Cap the panel and follow the normal panel login procedure.

Phoenix Inoculum Broth

Phoenix Inoculum Broth can be used to make the initial McFarland suspension of microorganisms when utilizing Phoenix MIC only panels (PMIC, NMIC, SMIC). The Phoenix Inoculum Broth is filled with 2.2 mL of Phoenix ID Broth and will reach the correct inoculum density by using fewer bacterial colonies.

- 1 Follow Steps 1 – 9 (**General Panel Preparation**) to prepare the suspension of bacteria using the Phoenix Inoculum Broth rather than the Phoenix ID Broth.
- 2 For NMIC and PMIC panels, add one drop of the Phoenix AST Indicator to the Phoenix AST broth tube.

For SMIC panels, add one drop of the Phoenix AST-S Indicator to the Phoenix AST-S broth tube.
- 3 Transfer 25 μ L (50 μ L if low inoculum option is used) of the Phoenix Inoculum Broth suspension to the Phoenix AST or AST-S tube(s).

- 4 Follow Steps 16 – 19 (General Panel Preparation). Follow the normal panel login procedure.

Phoenix MIC (NMIC, PMIC, SMIC) Panels

Inoculum can be prepared using either Phoenix ID Broth or Phoenix Inoculum Broth (Step 4, General Panel Preparation).

- 1 Follow Steps 1 – 9 (General Panel Preparation) to prepare the suspension of bacteria.
- 2 For NMIC and PMIC panels, add one drop of the Phoenix AST Indicator to the Phoenix AST broth tube.

For SMIC panels, add one drop of the Phoenix AST-S Indicator to the Phoenix AST-S broth tube.
- 3 Transfer 25 μ L (50 μ L if low inoculum option is used with NMIC or PMIC panels) of the Phoenix ID Broth to the Phoenix AST or AST-S tube. Cap and gently invert. Retain the ID tube or Inoculum tube for an optional purity check (see above).
- 4 Follow Steps 16 – 19 (General Panel Preparation). Follow the normal panel login procedure.

4.4.1 Quality Control

Refer to the panel package insert for information on ID/AST Quality Control. Note that QC panels cannot be marked as critical. Refer to Figure 4–2.

Quality Control testing is recommended for each lot of panels. The QC Lot Support feature can facilitate QC panel tracking and testing. If the QC Lot Support feature is enabled, then the panel lot number must be defined prior to logging in QC panels.

- 1 Inoculate a panel with one of the organisms listed in the package insert. All microbial cultures, including QC organisms, are potentially infectious and should be treated with universal precautions. QC organisms are prepared for panel inoculation as specified in Section 4.4 Preparing Panels, immediately preceding this section.



NOTE

For most reliable results, it is recommended that the QC organisms be subcultured at least twice on two consecutive days onto TSA with 5% Sheep Blood agar before use in the Phoenix system. Use only well-isolated colonies.

For Yeast ID QC, the use of Sabouraud Dextrose Agar as the subculture media is also acceptable.

- 2 Log the panel in as a QC panel as follows:

- a Press the "panel login" soft key.

- b Press the "QC panel" soft key.

- c In the **Sequence #** field, type in or scan the panel's sequence number.
- d If desired, in the **Accession #** field, type in or scan an accession number. Press  to advance to the next field.
- e If a combination panel is being used, both the ID and AST checkboxes are automatically checked (enabled).

14

PANEL LOGIN
GNIDAST

Sequence #

Accession #

Isolate #

Test Strain

☐ ID

☐ AST

Tech ID

Panel Lot #

ID Broth Lot #

AST Broth Lot #

Indicator Lot #

Exp. Date

Exp. Date

Exp. Date

Exp. Date

Figure 4-2 – QC Panel Login Display

For Yeast ID panels, the Media field appears to the right of the Isolate # field.

- f The **Isolate #** field defaults to isolate number 1. Type in the isolate number, or press UP ARROW or DOWN ARROW to increase or decrease the number. Valid isolate numbers are 1 to 20. Press **Tab** to advance to the next field.
- g For Yeast ID panels, you must specify a media type in the **Media** field. If you do not specify a media type, a workflow error is generated when you try to save the panel. If you do not login a Yeast ID panel before placing it in the instrument for testing, the panel aborts after the first reading because no media has been specified.

You can configure a default media type (see Section 2.3.1.2 – More Configuration Options), which appears when a Yeast ID panel sequence number is scanned during login.

To select a different media, press DOWN ARROW to drop down a box listing all media types (abbreviations), sorted alphabetically. (When you highlight the media type, the full name appears at the top right of the display.) Use UP ARROW or DOWN ARROW to highlight the desired media. Press the "select" soft key or **Enter** to select that media.

- h In the **Test Strain** field, press DOWN ARROW to drop down the QC organism selection box. (Only the predefined, required test strains appear [if available], otherwise all of the test strains appear. If a Sequence number has not been entered, "No Data Available" is displayed.) Use UP ARROW or DOWN ARROW to highlight the desired QC organism. Press **Tab** to select the organism. Press **Tab** to advance to the next field.
- i In the **Tech ID** field, enter the identification for the technologist performing the QC

test. Up to 3 alphanumeric characters are accepted. Press  to advance to the next field.

The following two fields do not appear when EpiCenter communications is enabled.

- j In the **Panel Lot #** field, type in or scan the panel's lot number. Lot numbers must be 7 numeric digits. Press  to advance to the next field. (If the QC Lot Support feature is enabled, this field is completed automatically when the Sequence # barcode is scanned.)
- k In the first **Panel Lot # Exp(iration) Date** field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date. (If the QC Lot Support feature is enabled, this field is completed automatically when the Sequence # barcode is scanned.)
- l If desired, in the **ID Broth Lot #** field, type in or scan the broth lot number. Lot numbers must be 7 numeric digits. Press  to advance to the next field.
- m If desired, in the first **ID Broth Lot # Exp(iration) Date** field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.
- n If desired, in the **AST Broth Lot #** field, type in or scan the broth lot number. Lot numbers must be 7 numeric digits. Press  to advance to the next field.
- o If desired, in the first **AST Broth Lot # Exp(iration) Date** field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.
- p If desired, in the **Indicator Lot #** field, type in or scan the broth lot number. Lot numbers must be 7 numeric digits. Press  to advance to the next field.
- q If desired, in the first **Indicator Lot # Exp(iration) Date** field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.
- r Press the "save" soft key to save the information.

- 3 Place the panel in the instrument (refer to Section 4.6).

- 4 When panel testing is complete, review the results for accuracy in the panel results display (see Section 4.10).

4.5 Panel Login

4.5.1 Using the Instrument Interface

The Liquid Crystal Display (LCD) presents all the information needed to monitor instrument status, to enter and remove panels, set up the instrument, print reports, and perform routine instrument maintenance. The information is presented in the form of screens and icons that graphically represent the information (e.g., a thermometer indicates the current temperature). The top region of the main status display presents instrument status information that is updated every few seconds. The middle region of the display initially presents station status information. As you access different screens and functions, they appear in place of the top and middle status windows. At the bottom of the displays, a series of icons shows the current software (soft) key definitions. Display regions are discussed in greater detail in **Section 5 – Reference**.

Most of the operations you perform at the instrument are initiated by pressing soft keys. Soft key definitions change as you select their related functions and access different displays and operations.

The keyboard can be used instead of the soft keys. The soft keys, from left to right, correspond to the keyboard function keys **F1** – **F8**. The UP ARROW, DOWN ARROW, and TAB keys correspond to the same keys on the keyboard.

In addition to the soft keys, six hard keys, themselves marked with screened icons, perform fixed functions regardless of the current display or operation. The LOAD PANELS key is used to prepare the instrument for loading panels for testing; UNLOAD PANELS key is used to remove panels from the instrument; SILENCE ALARM key quiets any audible alerts; the UP and DOWN ARROW keys are used to scroll through lists or to increase or decrease values; and the TAB key is used to move the cursor to the next field on a display. Pressing the SILENCE ALARM + UP ARROW keys together prints the current display.

4.5.2 Logging in Panels

Depending on the type of panel selected, different fields appear on the Panel Login display. (For information on logging in QC panels, refer to Section 4.4.1.)

The inoculum density of the panel is set in Configuration. There is no way during Panel Login to change the density setting. The only way to use a different inoculum density is by blackening well A-17 as described in **Section 4.4 Preparing Panels, ID Inoculum Density Flexibility**.

For Yeast ID panels (only), in order to insure optimal system performance, you must provide the correct media type via a drop-down box or by using the default media setting. The media type selection only applies to Yeast ID panels and is not displayed for other panel types.

Log the panel into the instrument as follows:

- 1 Press the "panel login" soft key.

- 2 In the **Sequence #** field, type in or scan the panel's sequence number.
- 3 In the **Accession #** field, type in or scan an accession number. Press **Tab** to advance to the next field (not required if barcode is scanned).
- 4 The **Isolate #** field defaults to isolate number 1. Type in the isolate number, or press UP ARROW or DOWN ARROW to increase or decrease the number. Valid isolate numbers are 1 to 20. Press **Tab** to advance to the next field. (You must enter an isolate number if you have entered an accession number.)
- 5 For Yeast ID panels, you must specify a media type in the **Media** field. If you do not specify a media type, a workflow error is generated when you try to save the panel. If you do not login a Yeast ID panel before placing it in the instrument for testing, the panel aborts after the first reading because no media has been specified.

You can configure a default media type (see Section 2.3.1.2 – More Configuration Options), which appears when a Yeast ID panel sequence number is scanned during login.

To select a different media, press ▼ to drop down a box listing all media types (abbreviations), sorted alphabetically. (When you highlight the media type, the full name appears at the top right of the display.) Use ▼ or ▲ to highlight the desired media. Press the "select" soft key or **Enter** to select that media.

- 6 If you want to receive special notification (audible alarm and/or automatic printing of a Lab Report) when panel results are obtained (ID only, partial, or complete), press the "select" soft key to check the **Critical** checkbox. More information about critical panels is provided in Section 2.3.5.
- 7 If you are only using either the ID or AST portion of a combination panel, press the "tab" key to advance the cursor to the ID or AST checkbox. Press the "select" soft key to UN-check the part of the panel you are not using.
- 8 If you UN-check the ID checkbox, or if you are using an AST-only or Phoenix Plus panel, the **Organism ID** field appears. If the system is not performing the organism identification, you must provide the organism ID for SIR interpretation. (If an AST panel is being tested and you do not enter an organism ID, the panel will go to "Needs Attention" when the instrument completes reading. You will have to provide an organism ID in order for the BDXpert system to interpret MIC results.)

Press DOWN ARROW to select the desired organism from the drop-down box. Organisms are listed in alphabetical order. You can enter the first few characters of the organism name to jump to that portion of the list quickly. Use UP ARROW or DOWN ARROW to highlight the desired organism. Press the "select" soft key or **Enter** to select the organism. You can also scan the desired organism from the barcode list of organisms found in the Quick Reference Guide.

- 9 Press the "save" soft key to save the information.

- 10 Place the panel in the instrument (refer to Section 4.6).


The following soft key functions can be performed from Panel Login:

"save" – saves the information displayed

"clear" – clears the record currently displayed from the screen

"select" – checks/unchecks a checkbox field or selects highlighted item in drop-down list box

"QC panel" – brings up fields needed to login a QC panel

"repeat data" – enters the last accession number and media type for panels, or media type and lot information for QC panels as follows: if QC Lot Support is disabled (Panel Lot plus Expiration, ID Broth Lot plus Expiration, AST Broth Lot plus Expiration, Indicator Lot plus Expiration), if QC Lot Support is enabled (ID Broth Lot plus Expiration, AST Broth Lot plus Expiration, Indicator Lot plus Expiration)

"LIS Orders" – brings up a window showing any order records that have been received from the LIS system (see Section 4.5.3 for more information)

4.5.3 Associating and Removing LIS Orders

The Phoenix instrument can save up to 200 Order records sent from an attached LIS system. These Order records can be accessed through a Panel Login – LIS Orders display. The LIS Orders display is accessible from the main Panel Login display via a soft key.

To access the LIS Orders display:

- 1 From the Main Status screen, press the “panel login” soft key.

- 2 Type in or scan the panel sequence number. (When a sequence number is entered, the system displays only LIS orders containing that panel type.)
- 3 Uncheck the ID or AST checkbox to disable that portion of a Combination panel, if desired.
- 4 From the Panel Login screen, press the “LIS Orders” soft key.

- 5 The LIS Orders display appears (Figure 4–3).
- 6 In the LIS Orders field, press ▼ to drop down a box listing LIS order records for that panel type, sorted alphabetically by accession number. The box lists the following information: accession number, isolate number, and organism or test strain (if available). Note that critical panels are shown in red text in the dropdown box.
- 7 Use ▼ or ▲ to highlight the Order record you want to associate to the panel. Press the “select” soft key or **Enter** to select that Order.
- 8 Press the “save” soft key to save the association. The LIS order is then removed from the LIS Order list.

Removing LIS Orders

To remove one or more LIS Order records from the LIS Order list, you can access the Remove LIS Orders window. This window lists all the LIS Orders by panel type. Order records are first marked, then deleted via soft key functions.

To remove LIS order records:

From the LIS Orders window, press the “more LIS functions” soft key


4.6 Inserting Panels in the Instrument

Panels are inserted with the reaction (and panel sequence barcode label) side facing the interior of the instrument. Refer to Figure 4–4.

WARNING

THE INSTRUMENT DOOR IS ELECTROMECHANICALLY LATCHED AND IS CONTROLLED BY THE INSTRUMENT SOFTWARE. NEVER ATTEMPT TO DEFEAT THE DOOR LATCHING MECHANISM, OR TO OPEN THE DOOR WHEN THE "UNLOCKED" ICON IS NOT DISPLAYED. SERIOUS INJURY CAN BE CAUSED BY THE ROTATING CAROUSEL. IF THE CAROUSEL IS NOT COMPLETELY STOPPED WHEN THE DOOR IS OPENED, IMMEDIATELY CONTACT BD FOR SERVICE. NEVER ATTEMPT TO ROTATE THE CAROUSEL MANUALLY OR SERIOUS INJURY MAY RESULT.

- 1 Press the LOAD PANELS key.

- 2 When the audible signal sounds (and "door unlocked" icon appears), open the instrument door.

- 3 Select a panel holder where there is no panel in place and no LEDs are illuminated. Place the bottom part of the panel in the panel holder.
- 4 Press downward.
- 5 Pivot the top of the panel back into the panel holder.

NOTE

Do not snap the panel back into the holder. This may result in splashing of the inoculum, which may cause inaccurate results.

- 6 Allow the panel to move upward into place.
- 7 Close the instrument door. If more panels need to be inserted than there are available holders in the current section, wait for a moment for the carousel to rotate to provide additional available holders, and repeat Steps 2 through 7.
- 8 The system performs an inventory scan to locate any newly inserted panels and reads the barcodes of these panels.


Figure 4-4 – Inserting Panels

4.7 Panel Testing

Panels placed in the instrument are continuously incubated at 35° C. The instrument tests panels every 20 minutes: on the hour, at 20 minutes past the hour, and again at 40 minutes past the hour.

Before each test cycle the current time is read. If the test cycle completes successfully the time is saved. A successful test cycle occurs when there are no carousel errors and the user does not pre-empt the test by initiating a Load Panels, Unload Panels, or Locate Panel operation. When the user initiates these operations, the current time is compared to the last test cycle's start time. If more than 30 minutes have NOT elapsed, the requested operation (Load Panels, Unload Panels, etc.) is performed. If more than 30 minutes HAVE elapsed since the last test cycle's start time, a test cycle must successfully complete before the user is permitted to perform a Load Panels, Unload Panels, or Locate Panel operation. The system alert/door unlocked LED is also off (if a system alert is not active) and the "cannot get into instrument" tone sounds when more than 30 minutes has elapsed since the last test cycle when a panel operation is requested.

After each test:

- The summary counters on the Main Status Screen are adjusted to indicate current statuses
- Panels/records that require user action have a Needs Attention flag set
- System alerts are reported in the System Alerts list
- Auto Association occurs (see Section 4.9).


The LCD screen displays an icon to inform you when a panel test is taking place.

The number next to the icon indicates the number of minutes remaining in the test cycle.

4.8 Responding to Needs Attention Panels

As system operation progresses, some panels may enter a state where they are flagged as needing attention. Usually, you can perform some action to correct the situation. In some cases, you are informed that the condition warrants carefully reviewing the panel's results. A listing of "Needs Attention" panels can be accessed by pressing the "needs attention" soft key on the Main Status Screen.


The Needs Attention conditions (refer to Section 5.11) are, in order of priority:

- | | |
|--|---|
| • Test Aborted | Active soft keys: "ignore" "locate panel"
"delete panel" |
| • Can Not Identify Barcode | Active soft keys: "ignore" "locate panel" |
| • Can Not Read Panel Wells | Active soft keys: "locate panel" "delete panel" |
| • Panel Lot Expired | Active soft keys: "ignore" "clear NA" |
| • Invalid AST Results | Active soft keys: "ignore" |
| • Panel Missing | Active soft keys: "ignore" "delete panel" |
| • No Growth on Panel | Active soft keys: "ignore" "locate panel"
"delete panel" |
| • Panel Lot Undefined | Active soft keys: none to address condition |
| • Review QC results | Active soft keys: "panel results" "ignore"
"locate panel" "delete panel" |
| • Missing Accession Number | Active soft keys: "panel results"
"locate panel" |
| • Missing Organism ID | Active soft keys: "panel results" |
| • Can Not Determine Organism ID | Active soft keys: "panel results" "ignore" |
| • Invalid Organism ID | Active soft keys: "panel results" "ignore" |
| • Organism ID Conflict | Active soft keys: "panel results" |
| • BDXpert Rule Flagged (if enabled and manual) | Active soft keys: "panel results" |
| • Pending Too Long | Active soft keys: "ignore" |

The Needs Attention display lists the first 30 panels with a Needs Attention flag, sorted by reason code (in the same order as the list above) and then by accession number within each reason code. Panels without an accession number are listed first within each Needs Attention reason code. Note that even after panels have been removed from the Needs Attention list, you can still view the highest priority reason code in the Panel Results display.

Depending on the Needs Attention reason, the following soft keys (and associated keyboard key) can be used in some cases to resolve the condition.

	"panel results"		"ignore"
	"delete panel"		"locate panel"
	"print"		"toggle display"
	"clear NA"		

These soft keys are designed to enable you to correct (whenever possible) the condition causing the panel to need attention. Where it is not possible, other resolutions are provided via soft keys.

"Panel Results" takes you to the Panel Results display, where you may be able to correct the Needs Attention code by adding information to the panel record.

"Ignore" does not correct the Needs Attention reason code, but it informs the system that you are aware of the problem.

"Delete panel" takes you to Panel Results display where you can delete the panel.

"Locate panel" causes the carousel to rotate to the panels location and lights the station where the panel resides.

"Print" enables you to print a Needs Attention List report. The report is filtered if that view is selected, or unfiltered if that view is selected.

"Toggle display" switches the display between filtered on highest priority reason (the default view) and unfiltered display (showing all needs attention reasons listed separately). When the list is filtered, only the highest priority needs attention condition is shown in the list. Panels are shown in blue if additional needs attention conditions exist for that panel.

"Clear NA" clears the Needs Attention condition

To resolve panels that need attention:

- 1 Access the Needs Attention display.
- 2 Refer to Section 5.11 and the chart below for detailed information on the particular needs attention reason.
- 3 Press the appropriate soft key to access Panel Results (to add or modify information), delete the panel, ignore or clear the reason, or to find the panel in the instrument.

Condition	Meaning	Possible Cause(s)	Resolution(s)
Test aborted	A condition occurred which caused the panel to be invalid	- Ongoing panel not tested for more than 1 hour - Instrument turned off for more than 1 hour - Instrument door open more than 1 hour - Panel moved to a different tier/instrument - Incubator temperature too high or too low - System software did not execute testing algorithms for more than 1 hour - Media type not specified for Yeast ID panel	- Delete the panel - Repeat testing.
Cannot identify barcode	Internal barcode scanner could not read a panel barcode in a station where the instrument could determine that a panel was present	- Barcode label obscured or missing - Unknown panel type was placed into the instrument	Locate the panel in the instrument and examine the barcode - If the barcode is obscured, the panel must be discarded and another inoculated - If the barcode appears to be intact, replace the panel and close the door. After the next inventory, check the Panel Needs Attention screen. If the panel does not appear, the internal scanner can now read the panel
Cannot read panel wells	Internal barcode scanner has read a sequence number in a station but the instrument does not detect that a panel is present in that station	Panel not seated properly	Locate panel, remove it and replace it, as above
Panel lot expired	A panel was logged in (or has a test start date) with a panel lot number that has already expired	Panel being logged in or placed in the instrument is from an expired lot	Discard panel, reinoculate isolate using unexpired panel lot
Invalid AST results	At least one MIC cannot be interpreted. (Excludes QC panels.)	Refer to Table 5-1	Repeat testing of the antibiotic that cannot be interpreted
Panel missing	Internal barcode scanner read a sequence number on an Ongoing panel, but the panel is missing	- Panel removed before the test was completed - Internal scanner failure and the sequence number can no longer be read	If the panel is replaced on the same tier within 1 hour after removal, testing will resume. If the panel is not replaced, testing will abort
No growth on panel	No growth in growth control well. (Excludes ID only and QC panels.)	Instrument did not detect growth in the growth control well of the panel	- Subculture the organism (to insure that it is viable) and inoculate a new panel. Panel has been aborted - From the Needs Attention screen delete or ignore the panel
Panel Lot Undefined	A panel has been entered whose lot number is undefined (non QC panel)	Panel is from an undefined lot	Define panel lot or ignore
Review QC results	Status of a completed QC panel is "Review"	QC panel which yielded an incorrect ID or incorrect AST result for at least one antibiotic, or has no growth in growth control well	Repeat QC organisms. Check: culture purity, inoculum density
Missing accession number	Panel is missing accession or isolate information. (Orphan panel.)	Failure to enter accession or isolate information	Press the "panel results" key. Add the accession number using the barcode reader or by typing it in
Missing organism ID	Panel has no organism ID. (ID required to determine SIR results. Excludes QC panels.)	- For AST only panel, no ID has been entered - The panel has an unresolved tie or triplet instrument ID and has no related panel with an ID	Key in the organism ID (cursor will be at the ID field and the field will be highlighted). OR open the drop down box by pushing the DOWN ARROW key and keying in the first letters of the organism name. The list scrolls to those letters in the organism list. Press the "select" soft key to select the highlighted organism. Any BDxpert rules triggered by the given ID will automatically be presented at this point. Select or ignore the rules. When completed, push the "exit" key to see complete test results.

Table 4-2 – Needs Attention Resolutions

Cannot determine organism ID	Panel has an Instrument ID of "No identification" or has a related panel with "No identification" as a final ID	Panel has been in test for 12 hours and the instrument cannot determine the identification	Repeat testing. Check the following: • Culture purity • Inoculum density • Correct panel used? • Organism may not be in the Phoenix database
Invalid organism ID	Organism ID is not in Phoenix database	Panel received download information of an organism ID that is not in the Phoenix database	Use alternate method
Organism ID conflict	Panel completes testing and has at least one related un-finalized panel that contains a different Final ID. Excludes QC panels.	Completed panel has at least one related panel that contains a different ID.	Select "Panel Results" soft key and choose an organism. Selecting the organism may trigger BDXpert rules. If the rules are configured as manual, the "BDXpert" soft key appears.
BDXpert Rule triggered	Panel triggered at least one BDXpert rule and the rule is manually enabled in Configuration.	A BDXpert rule needs to be invoked in order to determine AST results. A panel is flagged if BDXpert rules are configured as enabled/manual. (Rules that are configured as enabled/automatic will automatically "trigger" and the panel will not be displayed in Needs Attention.	Each Expert rule is displayed individually, in sequence. To accept the rule, press the "OK" soft key; to reject the rule press the "reject" soft key. Use the "Refresh" soft key to delete BDXpert system "decisions" and start over. When all rules have been displayed and dispositioned, complete AST results show. After all rules have been invoked, press "Special Messages" soft key (if present) to view Special messages about characteristics of the organism.
Pending too long	Panel has not been scanned (during an inventory count) within 30 minutes of logging in Panel Login.	Panel was logged into the instrument but was not placed in the instrument within 2 reading cycles (approximately 30 minutes).	Repeat testing. Delete the panel.

Table 4-2 – Needs Attention Resolutions

4.9 Automatic Association of Panels

4.9.1 Overview

Automatic Association, or Auto Association, is a process that the system initiates after certain trigger events occur (described later). Its purpose is to associate ID results from one panel (an ID or Combination panel) to the AST panels that are related to it and that lack an ID. Panels are related by virtue of having the same accession and isolate number. (Thus panels with the same accession number and different isolate numbers are NOT related, and Auto Association does not occur between such panels.)

Auto Association often involves the initial stage of Organism ID Conflict Checking. The circumstances where Conflict Checking occurs are described below. Conflict Checking is the process where the system verifies that the panel which has triggered the process does not have a related panel with a Final ID that is different from its own. If there is a conflict, Auto Association does not occur between related panels.

Auto Association can help eliminate unnecessary work for you, eliminating the need to manually enter an organism ID for related panels. In addition, the Conflict Checking function helps ensure that the same organism ID is made for all the panels related to a patient specimen.

Note that Auto Association does not set the critical panel attribute for related panels.

4.9.2 Panel Types

The panel types whose ID information is associated to related panels are: ID panels, Combination panels (with both ID and AST sides enabled), and Combination panels with only the ID side enabled. Auto Association uses the Final ID.

The related panels whose IDs are set by Auto Association are: AST panels, and Combination panels with only the AST side enabled.

4.9.3 When Auto Association Is Not Performed

Auto Association is NOT performed in the following circumstances:

- ◆ When an organism ID conflict exists within an accession/isolate
- ◆ To a related AST panel that has been finalized or deleted, or that already has a Final ID
- ◆ Between related ID-type panels
- ◆ Between related AST-type panels
- ◆ Between QC and regular panels
- ◆ Between related QC panels
- ◆ If the Instrument ID field contains a Tie/Triplet, this information is not associated to any related AST panels, however the Final ID information IS auto associated

4.9.4 Organism ID Conflict Checking

Organism ID Conflict checking is performed using the Final ID for panels within the same accession/isolate (i.e., related panels) that are not finalized (and are not QC panels). If both panels have a Final ID (both are not blank) and they do not match exactly, an Organism ID Conflict is considered to exist. When an Organism ID Conflict exists, no Auto Association is performed. All the panels involved in the conflict have a Needs Attention flag set indicating that there is an organism conflict. Once you resolve the conflict (and provided that all other conditions for Auto Association are met), then Auto Association does occur.

Note that Organism ID Conflict checking is not performed on QC panels.

4.9.5 Auto Association and Related Actions

The trigger events that invoke Auto Association include: an ID or Combination panel status transitioning to Removable or accession/isolate information being modified. If there are no unresolved tie/triplet conditions and no Organism ID Conflicts exist:

- 1 The instrument searches for the first *unfinalized* related ID or Combination panel with a Final ID. If a panel with these criteria is found, the instrument performs Auto Association.

If a panel with the above criteria is not found, the instrument then searches for the first *finalized* related ID or Combination panel. If a finalized panel is found, the instrument performs Auto Association.

If there are related ID panels with no Final ID when an AST panel invokes Auto Association (and no Organism ID Conflict condition exists), the instrument sets a Missing Organism ID Needs Attention reason for the AST panel if its status is Removable.
- 2 Auto Association is performed for any related AST panels that do not have information in their Final ID fields if the AST panel's status is either ongoing or Removable (Auto Association is not performed on AST panels with a status of Pending).

Auto Association is the action of modifying an AST panel's Final ID field (when this field is blank) to contain the same information as a related ID panel's information.

The instrument ID information is only associated provided it is a single Organism ID and not a Tie/Triplet condition.

A Combination panel with both sides enabled uses its ID result to perform interpretations on the AST side of the panel.

- 3 When Auto Association is performed successfully, any existing Missing Organism ID Needs Attention reasons are automatically cleared from the AST panels.

4.9.6 Typical Auto Association Examples

Instrument ID Is Associated

ID and AST panels are logged with the same Accession/Isolate number. (The AST panel is not logged in with an Organism ID.) They are both placed into the instrument at the same time. The ID panel finishes testing. The panels look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A
AST		

When the ID panel has completed testing, it looks for related panels. It sees the AST panel with no Final ID. The ID panel checks for Organism ID Conflicts. There is no Organism ID Conflict condition, so the ID panel associates its Instrument ID to the AST panel record's Instrument ID field (regardless of the AST panel's status). The panels look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A
AST	Org A	Org A

Instrument ID and Final ID Are Associated

The ID panel is logged in with an Accession/Isolate number. It is placed into the instrument to test. It completes testing and looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A

The user changes the Final ID for this panel. It now looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org B

The user then logs in an AST panel with the same Accession/Isolate number as the ID panel. (The AST panel is not logged in with an Organism ID.) The AST panel completes testing and looks for related panels. It sees the ID panel and checks to make sure that there is no Organism ID Conflict within the Isolate. There is no Organism ID Conflict, so the AST panel associates the ID panel's ID to itself. The panels now look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org B
AST	Org A	Org B

Final ID Is Associated

An ID panel is logged in with an Accession/Isolate number. The panel is placed in the instrument to test. The ID panel finishes testing first with a Tie/Triplet condition. When this situation occurs, the panel looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	

When the Tie/Triplet is resolved (by the user selecting one of the organisms listed or another organism), the ID panel will look like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	Org A

The user then logs in an AST panel with the same Accession/Isolate number as the ID panel. (The AST panel is not logged in with an Organism ID.) The AST panel completes testing and looks for related panels. It sees the ID panel and checks to make sure that there is no Organism ID Conflict within the Isolate. There is no Organism ID Conflict, so the AST panel associates the ID panel's Final ID to itself. The panels now look like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	Org A
AST		Org A

4.10 Obtaining Results

The primary means of viewing results of panel testing is the Panel Results display. Panel Results enables you to perform the following functions:

- display a panel whose data is stored in the Phoenix database
- modify the information for a panel that is resident in the database
- mark a panel as critical
- print a Lab report on any panel that can be recalled or displayed on the screen
- locate a panel resident in the same instrument where the panel data is recalled
- delete panel information from the Phoenix database
- answer or display any triggered BDXpert Rules
- display any special messages
- finalize a panel

You can access the Panel Results display in several different ways:

- scanning a known panel sequence number or accession number with the external scanner while the Main screen is displayed
- selecting the "panel results" soft key (or its keyboard equivalent) from the Main Status Screen

- scanning or entering a known panel sequence number while the Panel Login screen is displayed
- selecting the "delete panel" soft key (or its keyboard equivalent) or the "panel results" soft key (or its keyboard equivalent) for a panel in the Needs Attention screen when the Delete option is active
- selecting the "panel results" soft key (or its keyboard equivalent) from the Batch Finalization screen
- selecting the "panel results" soft key (or its keyboard equivalent) from the Panel Inventory or Panel Lot Definition displays

Typical Panel Results displays are shown in Figures 4–5 and 4–6.

When a panel is recalled in the Panel Results screen, the following fields are displayed (note that all fields may not be displayed depending on the type of panel and whether the instrument is connected to an EpiCenter system):

- Accession Number
- Sequence Number
- Isolate Number
- Media
- Status
- Final ID (Test Strain for QC Panel)
- Panel Lot # (if QC Lot Support feature is enabled)
- (Panel Lot) Expiration Date (if QC Panel, and QC Lot Support feature is enabled)
- Test Start
- Test End
- Location
- Finalized Status (QC Status for QC Panel)
- Critical Status
- Needs attention flag
- Resistance marker flag
- Instrument Organism ID which may consist of:
 - Instrument Organism ID
 - Confidence Value
 - Up to 5 Supplemental Tests
- Biochemical Results (Figure 4–5) which consist of:
 - Biochemical Abbreviation
 - Actual Result (+, –, ?, or X)
 - Expected Result (+, –, V, or blank)
- AST Results (Figure 4–6) which consist of:
 - Antimicrob(ial)
 - Minimum Inhibitory Concentration (MIC)
 - Instrument SIR (I)
 - BDxpert SIR (E)
 - Final SIR Result (F)

And for QC panels:

- Tech ID
- Panel Lot # and Exp. Date
- ID Broth Lot # and Exp. Date
- AST Broth Lot # and Exp. Date
- Indicator Lot # and Exp. Date

0.5 **PANEL RESULTS**
GNIDAST 951

Accession #	1313	Test Start	10/21/03 14:27
Sequence #	429510001405	Test End	10/21/03 14:27
Isolate #	1	Location	1/A01
Status	COMPLETE	Modify Related Panels	☒
Final ID	Prot. mirabilis	Finalized	☒
		Critical	☒

Instrument	Organism ID	Confidence
	Prot. mirabilis	97%

Biochemical	Actual	Expected
A_ARARR	+	+
A_GLPRB	-	+
A_GLYB	+	+
A_GUGAH	-	+
A_LARGH	+	+

Icons: Save, Clear, Select, Media, Resistance Marker, Print, Exit

Figure 4-5 – Panel Results Display – ID Results

For Yeast ID panels, the Media field appears to the right of the Isolate # field.

The following icons signify the conditions noted:

A Needs Attention reason exists

A Resistance Marker has triggered

The following soft key functions can be performed from Panel Results:

"save" – saves the information displayed

"clear" – clears the record currently displayed from the screen

"select" – checks/unchecks a checkbox field or selects highlighted item in drop-down list box

0.5

PANEL RESULTS
NQC01

Accession #
Sequence #
Isolate #
Status
Final ID

Test Start
Test End
Location

Finalized ☐
Critical ☐

Antimicr MIC			Antimicr MIC			Antimicr MIC		
	I	E F		I	E F		I	E F
GM 4	S	S	ATM 4	S	R	CIP <=0.5	S	S
NN 8	I	I	AM >16	R	R	NOR <=0.5	S	S
IPM <=0.5	S	S	PIP >64	R	R	TE >16	R	R
CF >32	R	R	AMC 8/4	S	R			
FOX 2	S	S	TZP <=2/4	S	S			
CTX <=4	S	R	SXT <=0.5/9.5	S	S			
FEP <=2	S	R	FM <=32	S	S			

Figure 4-6 – Panel Results Display – AST Results

"find" – recalls the data for the accession # entered (or the panel results list display if there is more than one panel for the accession)

"toggle display" – switches between ID results, AST results, and Resistance Marker displays (if applicable) for combination panels

"toggle soft keys" – switches between the currently displayed set of soft keys and additional ones

"print Lab Report" – prints a copy of the Lab Report for the current panel

"locate panel" – causes the instrument to locate and indicate the current panel and unlock the door

"delete panel" – removes the current panel from the database

"go to BDExpert Rules" – accesses the BDExpert Rules Triggered display

"more antimicrobials" – displays additional antimicrobials if there are more than can be displayed on one screen

"special messages" – accesses special messages display

"exit" – exits the current display

4.10.1 Adding/Modifying ID Results

The system enables you to enter an organism ID manually, or to override the organism identified by the instrument. In addition, in some cases, the system will not be able to make a single identification determination based on panel results. In these cases, two or three organisms may appear in the Instrument ID field. When more than one organism appears as the Instrument ID, the system does NOT automatically enter an identification in the Final ID field. The user must select the desired organism. The actual organism may be determined either through supplemental tests, which are recommended in the Instrument ID window and/or the Special Messages display (accessible via the "special messages" soft key, shown above), or through other tests you perform.

To add/modify the Final ID:

- 1 From the Panel Results display, press the "tab" soft key to advance to the **Final ID** field.
- 2 Press the DOWN ARROW key to drop down the listing of organisms.
- 3 Type the first few letters of the organism name to jump to that alphabetical part of the organism listing.
- 4 Use the UP ARROW or DOWN ARROW to highlight the desired organism.
- 5 Press the  or **Spacebar** key, or the "select" soft key to select the highlighted organism.

0.25 **PANEL RESULTS**
PMIC/ID-100

Accession # Missing Organism ID
 Sequence # Test Start
 Isolate # Test End
 Status Location
 Final ID Finalized ☐
 Critical ☐

Instrument	Organism ID	Confidence	NVR	PXR	HGP	ANR
1	Staph. epidermidis	65%	0	100	bb	95
2	Staph. hominis	21%	0	0	-	5
3	Staph. saprophyticus	10%	100	0	-	95

Biochemical	Actual	Expected		
		1	2	3
A_ARARR	-	+	+	+
A_GLPRB	-	+	+	+
A_LALT	-	+	+	+

Figure 4-7 – Panel Results Display – ID Results Showing Triplets

4.10.2 Modifying AST Results

There may be times when you wish to modify the Final SIR results for a panel. Note that you cannot modify any Final SIR results if there are manual BDxpert rules pending. You must first accept or reject any pending rules, which allows the system to perform its final results processing. After the final processing is complete, you will be able to modify SIR results manually if desired.

To modify the Final SIR results:

- 1 From the Panel Results display, press the "tab" soft key to advance to the Final SIR field for the desired antimicrobial.
- 2 Press the DOWN ARROW to toggle the field entry among the following results:

S(usceptible)

I(ntermediate)

R(esistant)

X = Invalid, cannot interpret (see Table 5-1)

N(ot Susceptible)*

Blank (indicates ID is required)


* N indicates that the antimicrobial/organism does not have an upper breakpoint. When there is no upper breakpoint there are no criteria for calling an organism intermediate or resistant. This often occurs when there are no known resistant strains of an organism. In this case, if the MIC is below the lower breakpoint the SIR results can be reported as susceptible but if the MIC is above the lower breakpoint the only result that can be reported is N or not susceptible.

4.10.3 BDXpert Triggered Rules Display

By default, the BDXpert Rules are enabled and set to trigger automatically when test results are processed. However, you have the option in the BDXpert Rule Configuration to disable some rules (which means they would not be applied to results) or to set BDXpert rules to trigger manually. Any rules set to trigger manually must be accepted or rejected by the user before any additional rule processing or Final SIR determination occurs.

To access BDXpert Triggered Rules:

- 1 From the Panel Results display, press the "go to BDXpert rules" soft key to access the BDXpert Triggered Rules display.
- 2 The first pending rule is highlighted in the Rules (bottom left) window.
- 3 You can view the effect of the rule on Final SIR values in the AST Results (bottom right) window.
- 4 Press: the "accept" soft key to accept the pending rule; "cancel" to reject the rule; "re-execute" to cause all interpretation rules to be reapplied to the raw data; or "exit" to return to the Panel Results display.

4.10.4 Resistance Marker Display

If Resistance Markers have triggered on a particular panel, the "toggle display" soft key cycles from the AST results display to the Resistance Marker review display (Figure 4-8) to the ID results display.

The Resistance Marker review display shows the same information in the top window as the other two displays (ID and AST results). In the bottom window, the following information is shown:

- The BDXpert rule number that triggered the Resistance Marker
- The Resistance Marker code (abbreviation)
- The Resistance Marker Name
- The BDXpert rule Description

0.5

PANEL RESULTS

CT01P

Accession #	1313	Test Start	10/21/03 14:27
Sequence #	424940000111	Test End	10/21/03 14:27
Isolate #	1	Location	1/A03
Status	COMPLETE		
Final ID	Enteroc. faecalis		Finalized ☐ Critical ☐

Rule #	RM	Name
1507	VRE	Vancomycin Resistant Enterococcus
0224	HLGR	High Level Gentamicin Resistant
0227	HLSR	High Level Streptomycin Resistant

Rule Description

For isolates with an interpretation of resistant for high levels of streptomycin (HLSR), synergy with cell wall active agents (e.g., ampicillin, penicillin, vancomycin) is not achievable.

Figure 4-8 – Resistance Marker Review Display

4.11 Panel Finalization

The panel finalization function enables you to finalize panel test results singly or in a batch mode. You can finalize panel results for all panels that are ready for finalization or for individual panels. You can review individual panel results prior to finalization.

All fields on the Finalization display are read-only. To *modify* information for a panel, press the "panel results" soft key.

When the instrument is connected to the EpiCenter Data Management System, the Finalization display is not available.

To finalize/batch finalize panels:

- 1 From the Main Status screen, press the "finalization" soft key. The Finalization display (Figure 4–9) appears.

if there are no panels to be finalized, only the "exit" soft key is available.

If there ARE panels to finalize, the following soft keys appear:

"print summary report"

"finalize single"

"finalize all"

- 2 To finalize ALL eligible panels, press the "finalize all" soft key.
- 3 To finalize panels one at a time, press the "finalize single" soft key.

The first panel eligible for finalization is shown in the main body of the display. Press one of the following keys:


Figure 4-9 – Finalization Display

OK	"accept"	To finalize the current record and display the next record
	"cancel"	To skip the current record and display the next record
	"toggle display"	To show AST results (if ID results are displayed) or vice versa
	"panel results"	To go to the Panel Results display

- To print a summary report of all panels eligible for finalization, press the "finalization summary report" soft key. A sample is shown in Figure 4-10. This report shows the Accession # (primary sort), Isolate # (secondary sort), Test End date and time, Sequence #, Instr #/Station, and Finalized status (* if finalized, blank if not).

Finalization Summary Report					
08/15/01 09:11					
Accession #	Isolate #	Test End	Sequence #	Instr # / Station	Finalized
1313	1	08/15/01 13:13	420011234567	1/A02	*
1344	1	08/15/01 13:13	420027654321	1/B03	*
1489	1	08/15/01 13:13	420031234123	1/C04	*
1566	1	08/15/01 13:13	420047654765	1/D05	
End of Report					

Figure 4-10 – Finalization Summary Report

- 5 If you access the Panel Results display to add or modify information before finalizing a panel, be sure to save the modifications. Press the “exit” soft key (from Panel Results) to return to the Finalize Single display.
- 6 Once an individual (or batch of) panel(s) has been finalized, the “print Lab Report” soft key appears on the Finalization display. This key enables you to print Lab Reports for all panels that have been finalized during this session (up to 200 maximum).

- 7 Continue to review panel records and finalize until no additional panel records are displayed.

4.12 Responding to Alarms and Errors

As you perform activities at the Phoenix instrument, and as testing progresses, system alerts and errors may occur. Different types of alerts and errors are flagged by one or more of the following: "E" and "W" error codes, audible tones, the System Alert icon appearing on the LCD display, or the instrument's System Alert indicator flashing. Generally, the more serious the condition, the more ways the system notifies you of the problem.

CAUTION

When the system notifies you of alerts and errors, you should immediately respond to the condition.

System alerts, which comprise all "E" type error codes, are reported in the system alert list. These errors cause: the system alert icon to appear on the Main Status screen; the System Alert indicator to flash; and the audible alarm to sound. They can be reviewed by pressing the "system alert" soft key (see below) and scrolling through the listing with the UP/DOWN ARROW keys. The errors must be reviewed to clear the system alert condition.

Activity errors comprise all "W" type error codes. These errors (such as attempting to enter an invalid Sequence # in the Panel Login screen) cause the activity error icon to appear on the currently displayed screen. They do not put the system into an alert condition. These errors can frequently be cleared by simply performing the activity correctly (such as entering a valid Sequence #).

All the "E" and "W" error codes are discussed in detail in **Section 7.2 – Error/Alert Messages**. The audible tones are discussed in **Section 3.7 – Audible Tones and Alarms**. The audible alarm can be silenced by pressing the SILENCE ALARM key, if the alarm is a "system alert" type of alarm.

4.13 Printing Reports

If you have a compatible printer connected to the instrument, the following reports can be selected for printing (see Figure 4-11):

- ◆ Completed Lab Report – contains information for all panels whose status became “Complete” during the time period (up to the past 48 hours) you select.
- ◆ Accession Lab Report – provides information for a specified accession number.
- ◆ Needs Attention List Report – lists all panels in the instrument’s database that have a Needs Attention flag set that have not been ignored.
- ◆ Resident Panel Report – lists the panels contained in stations 1 – 25 for each tier.
- ◆ QC Lab Report – lists all Test Strain Organisms that have completed testing and all biochemical and/or antimicrobial MIC results (for a specified Panel Lot #) that exist in the Phoenix databases when the report is requested.*
- ◆ Cumulative QC Report – provides information on quality control testing of all panel types. *
- ◆ Daily Instrument Report – lists the status of the instrument at the time the report is generated, along with a maintenance checklist.
- ◆ Interpretation Rule Set Report – lists the antimicrobial breakpoints of the Interpretation Rule Set defined as the default Rule Set in the Instrument Configuration screen.
- ◆ BDxpert Rule Set Database Report – lists each BDxpert rule number and the text describing the rule, whether each rule is enabled/disabled and whether each rule is set to trigger automatically or manually in the system.
- ◆ Organism ID Code List Report – displays all Organism Names, abbreviations for all Organism Names that exist in the Phoenix database.
- ◆ Antimicrobial Code Report – displays all antimicrobials, abbreviations for all antimicrobials, and the LIS codes that exist in the Phoenix database from all panel configurations
- ◆ Lab Report – contains all information for a panel sequence number that exists in the Phoenix database, including all information in the Panel Results screen, any special messages, BDxpert Rules that triggered, or Needs Attention Reasons if they exist. This report can only be printed from the Panel Results display or Finalization display.
- ◆ Finalization Summary Report – contains all the panels eligible for finalization at the time the report was requested, as well as finalization status. This report can only be printed from the Finalization display.
- ◆ Panel Inventory Lab Report – prints a Lab Report for all panels in the instrument. This report can only be printed from the Panel Inventory display.
- ◆ Custom Breakpoint Difference Report – prints a listing of differences between old breakpoints and new ones after a Phoenix Update Disk or Install/Upgrade operation. This report can only be printed from the Custom Interpretation Rule Set (Configuration) display.
- ◆ Current QC Panel Lot Report – contains information on the most recent QC test for each of the required strains for a panel lot, up to a maximum of 20 strains. This report can only be printed from the Panel Lot Definition display.
- ◆ Historical QC Panel Lot Report – contains information on all tests for a strain for the current instrument, up to 200 tests. This report can only be printed from the Panel Lot Definition display.
- ◆ Panel Lot Report – lists all the panel records for any panel lot number in the current instrument.
- ◆ Panel Lot Database Report – lists all the defined panel lots in the current instrument, and provides statistical and reference information on those lots.

*Report does not appear when EpiCenter communications is enabled.

More information on these reports and samples are provided in **Section 5 – Reference**.

PHOENIX REPORTS	
Completed Lab Report Accession Lab Report Needs Attention List Report Resident Panel Report QC Lab Report Cumulative QC Report Daily Instrument Report Interpretation Rule Set Report BDxpert Rule Set Database Report Organism ID Code List Report Antimicrobial Code Report	Start Date/Time 11 / 13 / 07 16 : 00 Prior Period 24 Hours QC ☐

Figure 4-11 – Reports Menu

To print a report:

- 1 From the Main Status Screen, press the "print reports" soft key.

- 2 The report menu appears.
- 3 Use the UP ARROW or DOWN ARROW to highlight the report you want to print.
- 4 Complete any additional fields (such as an accession number for the Accession Lab Report) and press the "print reports" soft key to print the report. (See Section 5 for additional field requirements.)

Several reports can also be printed from the displays that relate to them (e.g., Needs Attention List Report).

4.14 Unloading and Discarding Panels


When panel testing is completed, panels may be removed from the instrument and discarded.

WARNING

THE INSTRUMENT DOOR IS ELECTROMECHANICALLY LATCHED AND IS CONTROLLED BY THE INSTRUMENT SOFTWARE. NEVER ATTEMPT TO DEFEAT THE DOOR LATCHING MECHANISM, OR TO OPEN THE DOOR WHEN THE "UNLOCKED" ICON IS NOT DISPLAYED. SERIOUS INJURY CAN BE CAUSED BY THE ROTATING CAROUSEL. IF THE CAROUSEL IS NOT COMPLETELY STOPPED WHEN THE DOOR IS OPENED, IMMEDIATELY CONTACT BD FOR SERVICE. NEVER ATTEMPT TO ROTATE THE CAROUSEL MANUALLY OR SERIOUS INJURY MAY RESULT.

To remove panels:

- 1 Press the UNLOAD PANELS key.

- 2 When the "door unlocked" icon appears, open the instrument door.

- 3 All panels that are ready to be removed are indicated by a SOLID GREEN LED indicator.
- 4 Remove the panels by pushing the panel down, pivoting the top of the panel outward, and pulling the panel out of the panel holder.
- 5 If there are completed panels that are not in the accessible stations, close the door and allow the instrument to reposition the carousel to provide access to those panels. When the "door unlocked" icon appears (the carousel halted tone also sounds), open the door and continue removing completed panels.
- 6 Discard the panels in a biohazard container.

4.15 Installing New Panel Formats

From time to time, new panel formats will be designed that you may wish to use with your system. The Phoenix instrument is able to update its panel database via a series of barcodes, included with the new panels, that you scan. This allows you to update the database quickly and easily to use the new panels.

The printed barcodes provide the following information: the new panel's name; the system software version needed to use the new panel; the antimicrobials on the new panel; the concentrations of these drugs; and verification information to check that all the required data is entered, in the correct order, and is valid.

Only AST panels and the AST portion of Combination panels can be updated in this manner.

If you attempt to use a new panel format without first updating the instrument's panel database, an activity ("workflow") alert (W100 Invalid Panel Sequence Number) occurs, which prohibits you from logging in the new panel. Simply perform the panel update procedure described below to update the panel database, then you can login your new panels. The update process typically takes only a few minutes to perform.

To perform a panel update:

- 1 Display the Main Status screen.
- 2 Remove the BD Phoenix Panel Format page from the panel carton. Scan the first barcode on the page. The File Transfer display appears automatically.
- 3 You don't have to keep pressing the barcode trigger. Once the File Transfer operation begins, the beam remains on unless a scanning error occurs. Sweep the scanner down the first column (about 3 - 4 barcodes per second can be processed easily by the system).
- 4 When the first column of barcodes is complete, begin scanning the second column of barcodes. You **DO** have to scan the barcodes in order: first all the barcodes in the left column from top to bottom, then those in the second column, etc.
- 5 A graphical progress indicator on the display shows where you are in the update process. The current barcode number and name are also shown, for informational purposes.
- 6 **If any scanning errors occur**, a message appears on the display telling you which barcode number was not scanned. Simply press the scanner trigger and begin scanning from the numbered barcode onward. If any other type of error occurs, contact your BD representative for assistance.
- 7 To cancel the panel update operation, stop scanning the barcodes and press the "exit" soft key. The update is aborted and no changes are saved to the Phoenix instrument database.
- 8 When all the barcodes are scanned, the system processes the information and updates its panel database. Messages inform you of the progress of these operations. When the processing is complete, the message, "Merge complete" appears on the screen. Press the "exit" soft key to exit the File Transfer display.
- 9 When it's convenient, create a "panel data" disk with the new panel information you just added. See below, **Saving Your Panel Updates**.

Saving Your Panel Updates

After you have updated your panel database, we strongly recommend that you create a “panel data” disk which includes this new panel information. First, if there is a software update subsequently, the new panel database installed as part of *that* update may **not** contain the panel(s) you just added manually. A “panel data” disk enables you to easily “re-update” your panel database. Second, if you have multiple instruments, or a BD EpiCenter system, you can use the “panel data” disk to update *their* databases.

To create a “panel data” disk:

See Section 6.2.4.7 – Save Panel Configuration.

To install this “panel data” on a Phoenix instrument:

See Section 6.2.4.8 – Install Panel Configuration.

To install “panel data” files on a BD EpiCenter system:

Refer to the BD EpiCenter system help or user's manual.

4.16 Phoenix Update Disk (PUD)

Phoenix Update Disks may be provided from time to time to update Phoenix instrument databases and support files. These updates do not affect or change the basic instrument application software. For this reason, the Main Status Screen shows both the software version and the PUD version near the top of the display.

Instructions for updating the instrument using the PUD operation are provided in Section 6.2.4.9 – Phoenix Update Disk.

4.17 LIS Operations

4.17.1 General

The LIS communications feature enables the BD Phoenix instrument to exchange information with a compatible Laboratory Information System (LIS). LIS communications can be configured to exchange Order records and Results records at a variety of times.

Order records can be downloaded from the LIS system to the Phoenix instrument. These Order records can include the information listed in the second bullet below. If all this information is sent from the LIS system, then the panel is automatically logged in just as if the login were done at the instrument. If the panel sequence number is omitted, the Order record can be associated to a specific panel manually in Panel Login. Order records can be configured to be uploaded to the LIS system also.

Results records are uploaded from the Phoenix instrument to the LIS system. These records consist of:

- ◆ Header record (Delimiter fields, sender name, version number, message date/time)
- ◆ Order record (Accession number, Isolate number, Organism, Test ID, Sequence number, Priority, Report type)

Note that if an Organism code is included in the order record, for either a new or existing panel, the ID side of the panel is disabled.

- ◆ Comment record (contains Special Messages and/or BDxpert Rules)
- ◆ Results record
 - ◆ Panel sequence number
 - ◆ Instrument number
 - ◆ Instrument type
 - ◆ Instrument location (Station)
 - ◆ Time to result for identification or MIC produced
 - ◆ Test start time
 - ◆ Test end time
 - ◆ Test status (ongoing, complete, partial complete, complete with needs attention reason, complete with all ignored needs attention reasons, complete QC Pass, complete QC Review, pending, or rapid complete)
 - ◆ Result type
 - ◆ Antimicrobial code
 - ◆ MIC value
 - ◆ S / I / R / No Interp / Error value
 - ◆ Resistance marker
 - ◆ ID or Final ID
 - ◆ Results status (finalized or unfinalized)

The LIS Communications function is based on the American Society of Testing and Materials (ASTM) LIS Communications Standards (1381 and 1394), and is compatible with a number of popular LIS systems. For specific information on which LIS systems are compatible, contact your local BD representative. The *LIS Vendor Interface Specification* (available upon request) provides complete details on the BD Phoenix implementation of LIS Communications.

NOTE

If the Phoenix instrument is connected to a LIS, it cannot be connected to the EpiCenter™ system. However, if the Phoenix instrument is connected to EpiCenter, then a LIS connection can be established via EpiCenter.

4.17.2 Important Concepts

There are several concepts with which you should be familiar in order to understand the information presented here, and to insure that the LIS communications feature is set up properly for your laboratory. These concepts are discussed below.

LIS communication is able to send Results records from the instrument to the LIS (*upload*), Order records when panels are placed into the instrument, and queries to (and from) the LIS for Order records. Results upload can be configured to include or exclude Interpretation (SIR) results. The instrument can be set up to upload Results records only when the LIS requests them (*solicited upload*); or at one of the following *unsolicited upload* times: when panels are finalized; when panels complete testing or when complete panel records change; when ID or AST results are determined; or at a fixed time. QC panels and orphan panels are uploaded only when solicited by the LIS. If the instrument is configured for unsolicited upload, it still responds to requests from the LIS for results (solicitations or queries). If the instrument is configured for unsolicited uploads, the LIS must always be ready to receive data from the instrument.

By downloading barcode (panel) sequence numbers, the instrument can automatically login panels when records are received. If panel sequence numbers are omitted from the Order record, the panel can still be associated to the Order record in the Panel Login display.

If the panel sequence numbers are not included in downloaded records, the records that are downloaded are stored (up to 200 Order records) for manual association to panels. This process is described in **Section 4.5.3 – Associating and Removing LIS Orders**.

Organism Configuration and Antimicrobial Configuration enable you to enter the specific codes required by your LIS system for the organisms and antimicrobials uploaded in Results records. Refer to Sections 2.3.2.1 and 2.3.2.2 for additional information.

LIS configuration settings are independent of critical panel configuration settings. (E.g., if LIS configuration is set to send results only when the panel is complete, the results are **not** uploaded if the panel is critical and rapid reporting configuration is set for notification on ID results or partial results.)

4.17.3 Processing Panel Orders

Refer to **Section 4.5.3 – Associating and Removing LIS Orders**, for step-by-step instructions.

4.17.4 Routine System Operation

Operation of the Phoenix instrument with a LIS interface differs very little from routine operation of the system. The major difference is the ability to enter panel/accession data via the LIS into the system.

With LIS communications, patient information can be logged in at the LIS and transferred to the Phoenix instrument. (Consult your LIS manufacturer's operation manual for complete instructions on data entry and downloading records.)

LIS systems operate either in realtime mode, where the system automatically downloads each Order as it is logged in, or in batch mode where you log in and download multiple Orders in a group. After patient records are logged in, you should download them to the Phoenix instrument.

Any data sent to the Phoenix instrument that does not directly correlate to one of the fields defined as the "Order record" is ignored by the system. Any information sent from the LIS for a Finalized panel is rejected.

After the Order records have been downloaded, and the panels have been attached to those records (if necessary), routine system operation does not differ in any way. You should continue to load the instrument, print reports, monitor the system for complete panels, perform maintenance, etc. However, it is recommended that you be especially alert to and respond to any system or activity alerts that occur.

5 Reference

5.1 General

This section presents reference material on the Phoenix instrument user interface. All the screens, icons, reports, and functions in the user interface are described in the order in which they are accessed through the Main Status screen. The following information is presented:

- ◆ Main Status screen
- ◆ Panel Login
- ◆ Panel Results
- ◆ Finalization
- ◆ Configuration/Maintenance
- ◆ Reports
- ◆ Panel Inventory
- ◆ System Alerts
- ◆ Needs Attention

5.2 Software Tree

The following is a hierarchical list of all displays/functions in the system. The sections where these activities are discussed in detail are noted in parentheses.

- Main Status screen (Section 5.3)
 - Panel Login (Section 5.4)
 - Panel Results (Section 5.5)
 - BDXpert Rules Triggered (Section 5.5.1)
 - Resistance Markers (Section 5.5.2)
 - Special Messages (Section 5.5.3)
 - Lab Report (Section 5.5.4)
 - Finalization (Section 5.6)
 - Finalization Summary Report (Section 5.6)
 - Configuration/Maintenance (Section 5.7)
 - Instrument Configuration (Section 2.3.1)
 - General Instrument Configuration (Section 2.3.1.1)
 - Laboratory Information (Section 2.3.1.2)
 - Password Administration (Section 2.3.1.3)
 - Maintenance (Section 4.2, 2.5, 7.3, 6.2.2, 6.2.4)
 - Communications Configuration (Section 2.3.2)
 - Organism Configuration (Section 2.3.2.1)
 - Antimicrobial Configuration (Section 2.3.2.2)
 - Custom Interpretation Rule Set (Breakpoint) Configuration (Section 2.3.3)
 - BDXpert Rule Configuration (Section 2.3.4)
 - Rapid Reporting Configuration (Section 2.3.5)
 - Panel Lot Definition (Section 2.3.6)
 - Reports (Section 5.8)
 - Completed Lab Report (Section 5.8.1)
 - Accession Lab Report (Section 5.8.2)
 - Needs Attention List Report (Section 5.8.3)
 - Resident Panel Report (Section 5.8.4)
 - QC Lab Report (Section 5.8.5)
 - Cumulative QC Report (Section 5.8.6)
 - Daily Instrument Report (Section 5.8.7)
 - Interpretation Rule Set Report (Section 5.8.8)
 - BDXpert Rule Set Database Report (Section 5.8.9)
 - Organism ID Code List Report (Section 5.8.10)
 - Antimicrobial Code Report (Section 5.8.11)
 - Lab Report/QC Lab Report (Section 5.8.12)
 - Finalization Summary Report (Section 5.8.13)
 - Custom Breakpoint Difference Report (Section 5.8.14)
 - Panel Lot Report (Section 5.8.17)
 - Panel Lot Database Report (Section 5.8.18)
 - Panel Inventory (Section 5.9)
 - System Alerts (Section 5.10)
 - Needs Attention (Section 5.11)

5.3 Main Status Screen

The Main Status screen (Figure 5–1) is the default screen that is displayed when no activity has been initiated. This screen presents up to eight icons representing soft key definitions that allow you to perform activities. These icons are as follows:

- ◆ The “panel login” icon, where you can login new panels.
- ◆ The “panel results” icon, which allows you to view and modify the ID or AST results of panels or delete them.
- ◆ The “finalization” icon, where you can finalize panels that have completed testing. (Not available when instrument is connected to the EpiCenter Data Management System.)
- ◆ The “configuration” icon, where you can set up instrument, laboratory-specific rules, and communications parameters, as well as perform instrument maintenance checks.
- ◆ The “print reports” icon. This icon appears only if your printer is attached, turned on, and online.
- ◆ The “panel inventory” icon, which allows you to review a list of all panels in the instrument and sort the list by various criteria (such as by panel sequence number, accession number, result, etc.)
- ◆ The “system alert” icon, which allows you to review any system alerts that may have occurred or that may still exist. (Icon appears only if system alert is present.)
- ◆ The “needs attention” icon, which allows you to respond to panels that have been flagged as needing attention. Usually, you can perform some action to correct the situation. (Icon appears only if there are panels that need attention.)

The Main Status screen has three areas: the top region presents system status information; the middle region presents station status information; and the bottom region displays the current soft key definitions.

The system status information presented at the top of the screen includes the following:

Current date and time

Carousel temperature

Instrument state

idle

#

testing

##

inventory scan

#

processing test data

#

= number of minutes until next test

= number of minutes left in current test for testing state

adjusting UV light source

adjusting RGB light source

* = number of minutes adjustment has been in progress

When the “processing test data” icon is displayed, do not change any configuration settings.

Door unlocked (blank if locked)

Instrument number

Software version / Phoenix Update Disk (PUD) version (VX.XXY) / VA.AAB

Inoculum density (red for Gram negative [top], blue for Gram positive [bottom])

Network okay

Network down

LIS enabled

Station statuses presented in the middle area of the display include the following:

Ongoing

#

Removable

#

Available

#

Blocked

#

The numbers (#) to the right of the icons represent the number of stations/panels that exist for a status.

The soft key assignment area (very bottom of the screen) shows the current soft key definitions. These definitions change as you access different functions and screens. The Main Status screen icons include:

"panel login"		"panel results"	
"finalization"		"configuration"	
"print reports"		"panel inventory"	
"system alert"		"needs attention"	

Note that the "finalization" icon does not appear when the Phoenix instrument is attached and communicating with an EpiCenter system.


Figure 5-1 – Main Status Screen (with optional LIS indicator)

5.4 Panel Login

The Panel Login display enables you to log in demographic information about panels to be tested. Depending on the type of panel being used, not all fields listed below may appear on the display. More than one of the same panel type may be logged in for an accession number, but the system notifies you of the duplication with an activity error message.

The inoculum density of the panel is set in Configuration. There is no way during Panel Login to change the density setting. The only way to use a different inoculum density is by blackening well A-17 as described in **Section 4.4 Preparing Panels, ID Inoculum Density Flexibility**.

See Figure 5–2.

To access Panel Login:

Press the "panel login" soft key from the Main Status screen

or

Scan a new Accession with the handheld barcode scanner or scan a new panel barcode with either scanner from the Main Status screen

Panel Login soft keys:

"save"

"clear"

"select"

"QC panel"

"repeat data"

"LIS Orders"


"exit"

Panel Login fields:

Sequence

Type in or scan the panel's sequence number. (If you scanned the Sequence # to access Panel Login, this field is completed automatically.) The sequence number contains digits that identify the panel type. Valid sequence numbers are 12 digits long and must conform to BD panel sequence number specifications. You must enter a panel sequence number to save a record. If an existing sequence number is entered, the system takes you to the Panel Results display.

The screenshot shows a 'PANEL LOGIN' screen for 'NMIC/ID-1'. At the top left is a circular icon with a gear and the number '14'. The title 'PANEL LOGIN' is centered, with 'NMIC/ID-1' below it. On the right side, there are three checkboxes labeled 'Critical', 'ID', and 'AST'. The main area contains three input fields: 'Sequence #' with the value '42501777777', 'Accession #' with the value '1313', and 'Isolate #' with a dropdown menu showing '1'. Below these fields is a large empty rectangular area. At the bottom, there is a row of eight icons: a floppy disk, a hand pointing at a screen, a hand pointing at a button, a curved arrow, a bar chart, a 'LIS' label, and a door with an arrow pointing out.

Figure 5-2 – Panel Login Display

For Yeast ID panels, the Media field appears to the right of the Isolate # field.

Accession

Type in or scan an accession number, up to 20 alphanumeric characters. (If you scanned the accession barcode to access Panel Login, the Accession # is completed automatically.) Spaces at the beginning or end of the number are ignored, but spaces within the number are valid. The following characters CANNOT be used in accession numbers: * ? [] ! # | .

Isolate

Defaults to isolate number 1. Type in the isolate number, or press UP ARROW or DOWN ARROW to increase or decrease the number. Valid isolate numbers are 1 to 20. If you have entered an accession number, you must also enter an isolate number. If you change the isolate # from "1" to another value, you must enter an accession number.

Media

This field appears only if a Yeast ID panel sequence number is scanned or typed in.

Defaults to media type specified in Configuration (optional). To select a different media for the isolate, press ▼ to drop down a box listing all media types (abbreviations), sorted alphabetically. (When you highlight the media type, the full name appears at the top right of the display.) Use ▼ or ▲ to highlight the desired media. Press the "select" soft key or Enter to select that media. You can select or change a media type until the first panel test. The panel test is aborted if a media type is not specified before the first test cycle.

You must select a media type to save the panel information.

Critical Checkbox

Defaults to unchecked. Press the "select" soft key to check the box and mark the panel as critical. The system can be configured (Rapid Reporting configuration, Section 2.3.5) to notify you of critical panel results (ID results obtained, partial panel results, or complete panel results) by sounding an audible alarm and/or printing a Lab Report automatically (or neither).

ID Checkbox

Defaults to checked if a combination panel sequence number is scanned. If you are using only the AST portion of a combination panel, press the "select" soft key to UN-check the box. This checkbox appears only when a combination panel is used. This field is automatically checked (enabled) when you log in a QC panel.

AST Checkbox

Defaults to checked if a combination panel sequence number is scanned. If you are using only the ID portion of a combination panel, press the "select" soft key to UN-check the box. This checkbox appears only when a combination panel is used. This field is automatically checked (enabled) when you log in a QC panel.

Organism ID/Test Strain

Press DOWN ARROW to drop down the organism selection box. Organisms are listed in alphabetical order. You can enter the first few characters of the organism name to jump to that portion of the list quickly. Use UP ARROW or DOWN ARROW to highlight the desired organism. Press  to select the organism.

This field appears only for AST-only panels or when using the AST side of a Combination panel.

For QC panels, this field is named Test Strain, and lists only the strains of organisms predefined in the database, sorted by strain number (alphanumerically). A test strain must be entered to save a QC panel.

Tech ID

This field appears for QC panels. Enter the identification for the technologist performing the QC test. Up to 3 alphanumeric characters are accepted. A tech ID must be entered to save a QC panel.

Panel Lot

This field appears for QC panels. Type in or scan the panel's lot number. Lot numbers must be 7 digits. A lot number must be entered to save a QC panel. (When the QC Lot Support feature is enabled, this field is completed automatically when the Sequence # barcode is scanned.)

(Panel Lot #) Exp. Date

This field appears for QC panels. In the first Panel Lot # Exp(iration) Date field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date. An expiration date must be entered to save a QC panel. (When the QC Lot Support feature is enabled, this field is completed automatically when the Sequence # barcode is scanned.)

ID Broth Lot

This field appears for QC panels. Type in or scan the broth lot number. Lot numbers can be up to 7 characters.

(ID Broth Lot #) Exp. Date

This field appears for QC panels. In the first ID Broth Lot # Exp(iration) Date field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press **[Tab]** to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press **[Tab]** to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.

AST Broth Lot

This field appears for QC panels. Type in or scan the broth lot number. Lot numbers can be up to 7 characters.

(AST Broth Lot #) Exp. Date

This field appears for QC panels. In the first AST Broth Lot # Exp(iration) Date field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press **[Tab]** to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press **[Tab]** to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.

Indicator Lot

This field appears for QC panels. Type in or scan the indicator lot number. Lot numbers can be up to 7 characters.

(Indicator Lot #) Exp. Date

This field appears for QC panels. In the first Indicator Lot # Exp(iration) Date field, press UP ARROW or DOWN ARROW to enter the expiration day, month, or year (depending on your configuration). Press **[Tab]** to advance to the next date field and press UP ARROW or DOWN ARROW to enter the expiration day, month, or year. Press **[Tab]** to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date.

LIS Orders

This field appears when the "LIS Orders" soft key is pressed. Entering a Sequence number and pressing the DOWN ARROW opens a drop-down box listing up to 200 Order records received from the LIS. The list contains only those LIS order records for the same panel type as the Sequence number entered. The list is sorted alphabetically by accession number. It lists the following information: accession number, isolate number, and organism or test strain (if available). Any critical Order records are shown in red text in the box. Use the UP ARROW or DOWN ARROW to highlight the desired Order record and press the "select" soft key to associate that Order to the Sequence # entered.

5.5 Panel Results

The Panel Results display enables you to review and modify panel test results. The Panel Results display can be used for the following functions:

- ◆ display a panel whose data is stored in the Phoenix database
- ◆ modify the information for a panel in the Phoenix database

- ◆ mark a panel as critical
- ◆ print a Lab report on any panel that can be recalled or displayed on the screen
- ◆ locate a panel resident in the instrument
- ◆ delete panel information from the Phoenix database
- ◆ answer or display any triggered BDXpert Rules
- ◆ display any special messages
- ◆ finalize a panel

Panel Results (clinical panels) are retained for 31 days (possibly longer depending on number of QC panels tested). QC panel results are retained for at least 6 months.

Additional information about the panel is displayed in the display title area. The type of panel is shown below the title. Also, the inoculum density used for the identification (if applicable) is shown near the top left after the first test cycle completes.

If there is a Needs Attention flag for the panel, it appears above the **Test Start** field. If multiple Needs Attention reasons exist, only the highest priority reason is shown. If a Resistance Marker has triggered for the panel, an icon is shown just to the left of the **Location** field.

If a panel is recalled by accession number only, and there is more than one panel attached to the accession, the Panel Results List display appears. From this display, you can select the specific panel to review/modify. See Figure 5-3 for the Panel Results List display. See Figure 5-4a for a sample ID Panel Results display, and Figure 5-4b for a sample AST Panel Results display.

Note that depending on the type of panel being used, as well as other circumstances, not all fields listed below may appear on the display.

AST or Combination panel types (with at least the AST side enabled) each contain a set of antimicrobials. The instrument reports a result for each antimicrobial on the panel. Each antimicrobial reports an individual MIC value. Once an antimicrobial has a MIC value, the instrument calculates an Instrument susceptibility (SIR) value for each MIC value that has been determined. (The instrument requires the panel to have an Organism ID defined to interpret MIC values into Instrument SIR values. Also note that SIR values are not calculated by the instrument if it is attached to an EpiCenter Data Management System.)

If Rapid Completion is enabled, the instrument is able to provide BDXpert AST results (SIR) before determining actual MIC values. The instrument MIC values are provided as soon as they can be accurately determined. Within a test panel, some MIC values may be available earlier than others. The rapid completion feature can be used to predict resistance for uncompleted antibiotics using the ID alone (intrinsic resistance), or ID with completed MICs of related antibiotics, or resistance marker tests (BL, ESBL). The BDXpert system is used to make these predictions. This can be useful in situations where, for example, the results for drugs that have not yet received MICs would be of no clinical value based on the other results that are already available. Antimicrobials with Rapid Complete BDXpert interpretations are indicated by a "C" in the MIC column on Results screens and Lab reports.

When both a MIC value and Instrument SIR value have been determined for an antimicrobial, the instrument executes the BDXpert Rules (providing the BDXpert System is enabled). The instrument reports a value in the BDXpert SIR field if an enabled BDXpert Rule triggered and the reported BDXpert SIR value is different from the value in the Instrument SIR field.

To access Panel Results:

Press the "panel results" soft key from the Main Status screen

or

Scan an existing Accession with the handheld barcode scanner from the Main Status screen; or scan an existing panel barcode with either scanner from the Panel Login or Main Status screen

Panel Results can also be accessed in certain conditions from the Needs Attention and Batch Finalization displays.

Panel Results soft keys:

"save"

"clear"

"select"

"find"

"toggle display"

"toggle soft keys"

"print Lab Report"

"locate panel"

"delete panel"

"go to BDXpert Rules*"

"more antimicrobials"

"special messages"

"QC panel"

"exit"

* Note that this icon does not appear if the EpiCenter system is attached and communicating.

Panel Results fields:

NOTE

Modifying the Accession # or Isolate # may invoke Auto Association, which can change results.

Accession #

Type in or scan the accession number to recall. If only an accession number is entered, and there is more than one panel attached to the accession, the Panel Results List display appears (Figure 5-3).

This field can be modified for unfinalized panels. Enter up to 20 alphanumeric characters for the accession, excluding * ? [] ! # | . Modifying an accession number does not affect the accession number of any related panels.

Sequence #

Type in or scan the panel's sequence number to recall. This field cannot be edited.

PANEL RESULTS

Accession #

ACC333

Sequence #

Iso#	Panel Type	Final ID	Status	Test Start	Sequence #
1	GNIDAST	ENTBCLO	COMPLETE	02/15/00	420015555555
2	GNID		COMPLETE	02/15/00	420022222222
2	GNAST		COMPLETE	02/15/00	420031112774
3	GNIDAST	ESCCOL	COMPLETE	02/15/00	420016666666

Figure 5-3 – Panel Results List Display

A maximum of 30 panels can be displayed in the list. To select a panel, use the UP ARROW or DOWN ARROW to highlight the panel and press the "select" soft key.

Isolate

Press UP ARROW or DOWN ARROW to increase or decrease the number. Valid isolate numbers are 1 to 20. This field can be modified for unfinalized panels, however an existing isolate number cannot be changed to a blank number.

Media

This field is editable only before the first test cycle completes.

Shows the media type selected during Panel Login for Yeast ID panels only. The following values can appear in this field: blank for Unspecified Media Type; INVLD (Invalid Media Type); SABDX (Sabouraud Dextrose Agar); TSASB (Trypticase Soy Agar w/ 5% SB); COLSB (Columbia Agar w/5% SB); CHOC (Chocolate II Agar); SABEM (Sabouraud Dextrose Emmons); SABHI (Sab Brain Heart Inf Ag Deep)

Status

Read-only field shows the panel's testing status: Pending, Ongoing, Complete, or Rapid Complete (if enabled). Rapid Complete panels have not finished testing and show BDxpert results. If a Rapid Complete panel is removed from the instrument, it then becomes Complete. If it is left in the instrument, Rapid Complete MICs are replaced by actual MIC results as they are determined.

Final ID

Final ID can be completed automatically by the system from a single instrument-based ID or by the user selecting among tie instrument-based IDs. For ID panels (or ID-only portions of Combination panels) you cannot manually select an organism ID if there is no instrument-based ID listed.

Note that if the **Modify Related Panels** checkbox is checked, the ID in the database of all related panels (panels with the same accession and isolate number) is set with the same organism ID through the Auto Association function (see Section 4.9).

Press DOWN ARROW to drop down the organism selection box. Organisms are listed in alphabetical order. You can enter the first few characters of the organism name to jump to that portion of the list quickly. Use UP ARROW or DOWN ARROW to highlight the desired organism. Press  to select the organism.

When this field is modified by the user, the instrument changes the background color of the field to cyan. In addition, if the **Modify Related Panels** checkbox was checked at the time the **Final ID** was modified, all modified panels' **Final ID** fields are also displayed with a cyan background. When the field is modified, the system re-evaluates each antimicrobial's instrument SIR value, as well as re-evaluating BDxpert information.

The field cannot be modified if the panel is finalized.

For QC panels, this field is named **Test Strain**, and lists only the ATCC strains of organisms predefined in the database, sorted by strain number.

Panel Lot

Read-only field shows the panel lot number. This field shows only when the QC Lot Support feature is enabled.


  0.5		PANEL RESULTS GNIDAST 951																					
Accession #	1313	Test Start	10/21/03 14:27																				
Sequence #	429510001405	Test End	10/21/03 14:27																				
Isolate #	1 	Location	1/A01																				
Status	COMPLETE	Modify Related Panels	☒																				
Final ID	Prot. mirabilis 		Finalized ☒ Critical ☒																				
Instrument Organism ID		Confidence																					
Prot. mirabilis		97%																					
<table border="1" style="width: 100%; border-collapse: collapse;"> <thead> <tr> <th style="width: 30%;">Biochemical</th> <th style="width: 20%;">Actual</th> <th style="width: 20%;">Expected</th> <th style="width: 30%;"></th> </tr> </thead> <tbody> <tr> <td>A_ARARR</td> <td style="text-align: center;">+</td> <td style="text-align: center;">+</td> <td rowspan="5" style="text-align: center; vertical-align: middle;">  </td> </tr> <tr> <td>A_GLPRB</td> <td style="text-align: center;">-</td> <td style="text-align: center;">+</td> </tr> <tr> <td>A_GLYB</td> <td style="text-align: center;">+</td> <td style="text-align: center;">+</td> </tr> <tr> <td>A_GUGAH</td> <td style="text-align: center;">-</td> <td style="text-align: center;">+</td> </tr> <tr> <td>A_LARGH</td> <td style="text-align: center;">+</td> <td style="text-align: center;">+</td> </tr> </tbody> </table>				Biochemical	Actual	Expected		A_ARARR	+	+		A_GLPRB	-	+	A_GLYB	+	+	A_GUGAH	-	+	A_LARGH	+	+
Biochemical	Actual	Expected																					
A_ARARR	+	+																					
A_GLPRB	-	+																					
A_GLYB	+	+																					
A_GUGAH	-	+																					
A_LARGH	+	+																					
       																							

Figure 5-4a – Panel Results Display – ID

For Yeast ID panels, the Media field appears to the right of the Isolate # field.

Exp. Date

Read-only field shows the date and time that the panel lot expires. This field shows only for QC panels and only when the QC Lot Support feature is enabled.

Modify Related Panels Checkbox

This field appears when the Final ID field is modified. Defaults to checked when the Final ID is modified. If you do NOT want to modify the ID for unfinalized related panels when you modify the currently displayed panel, press the "select" soft key to UN-check the box. This checkbox has no affect on related QC panels. This checkbox does not appear if the EpiCenter system is attached and communicating.

Test Start

Read-only field shows the date and time that panel testing was started.

Test End

Read-only field shows the date and time that panel testing was completed.

Location

Read-only field shows the location of the panel, in the form I/Tst, where I is the Instrument number, T is the Tier, and st is the station number.

QC Status

This field shows the status of a QC panel. The field is blank until the panel status becomes Complete. Statuses are initially PASS or REVIEW. Review indicates that the panel has not passed. You should check any panels with Review status and determine why the panel did not pass.

The status is REVIEW if any of the following conditions occur:

- ◆ QC strain was identified incorrectly
- ◆ The test on a QC panel is aborted
- ◆ At least one of the antimicrobial results fail

From a status of REVIEW, you may set the final status to REPEAT if you determine that the panel failed due to error in preparation or handling of the panel. If it is not determined that a panel preparation/handling error was made, the final status should be set to FAIL. Selecting REPEAT or FAIL clears the "Review QC Results" Needs Attention condition.

Finalized

Press the "select" soft key or Spacebar to check the box and finalize the panel. The **Finalized** field is displayed as a read-only field when the Phoenix system is connected and communicating with the EpiCenter Data Management Center. In this case, all Finalization is done at the EpiCenter system. **Finalized** is not displayed for QC panels.

If you finalize a panel with Rapid Complete status, processing of MIC results stops, and the panel status becomes Complete. Any drugs that did not complete testing maintain their current MIC value (e.g., C in MIC column remains C, X remains X, etc.).

Critical

Press the "select" soft key or Spacebar to check the box and mark the panel as critical. The system can be configured (Rapid Reporting configuration, Section 2.3.5) to notify you of critical panel results (ID results obtained, partial panel results, or complete panel results) by sounding an audible alarm and/or printing a Lab Report automatically (or neither).

You can also unmark previously marked critical panels by unchecking the checkbox.

This field does not appear for QC panels.

Accession #		Test Start	
Sequence #		Test End	
Isolate #		Location	
Status		Finalized	
Final ID		Critical	
425190000030		10/10/00 10:00	
6		10/10/00 10:00	
COMPLETE		10/10/00 10:00	
Esch. coli			

Antimicr MIC	I	E	F	Antimicr MIC	I	E	F	Antimicr MIC	I	E	F
GM 4	S	S	S	ATM 4	S	R	R	CIP <=0.5	S	S	S
NN 8	I	I	I	AM >16	R	R	R	NOR <=0.5	S	S	S
IPM <=0.5	S	S	S	PIP >64	R	R	R	TE >16	R	R	R
CF >32	R	R	R	AMC 8/4	S	R	R				
FOX 2	S	S	S	TEP <=2/4	S	S	S				
CTX <=4	S	R	R	SXT <=0.5/9.5	S	S	S				
FEP <=2	S	R	R	FM <=32	S	S	S				

--	--	--	--	--	--	--	--	--

Figure 5-4b – Panel Results Display – AST

Antimicr(obial)

Read-only field showing the abbreviation for the antimicrobial.

MIC.

The Minimum Inhibitory Concentration value determined by the instrument. The following values may also appear in this field:

- | | |
|---|---|
| > | growth occurred for all concentrations of the antimicrobial |
| ≤ | no growth occurred in any of the concentrations of the antimicrobial |
| ? | MIC determination is pending (SIR values remain blank) |
| X | MIC value cannot be produced; or Final ID is not claimed in the Taxa listing (Appendix F.3) for AST testing; or the panel's drug dilution series does not cover the Phoenix reportable MIC range, or for additional causes, see Table 5-1. Check Special Messages display (Section 5.5.3) for an explanation. |
| C | Rapid Complete (MIC is pending, BDxpert SIR is based on ID and completed drug and/or resistance marker results). These values are replaced by actual MIC values as they are determined. |

For QC Panels, the following values may appear in this field:

Actual MIC value (number, ?, or X for error)

Expected MIC range (if defined)

QC Result

- P the actual MIC value is within the expected MIC range
- F the actual MIC value is not within the expected MIC range
- R Repeat: the actual MIC value is "X"
- blank no range is defined for the antimicrobial/organism combination

I(nstrument SIR)

This field is the instrument interpretation for the MIC based on the breakpoints currently running in the system (selected in Instrument Configuration).

The SIR value is blank for antimicrobials that require an ID to perform the SIR interpretation. In this case, the **BDXpert SIR** field (if present) and the **Final SIR** field are also blank.

The following represent the interpretation values:

- S Susceptible
- I Intermediate
- R Resistant
- N Not Susceptible*
- X Cannot produce interpretation

Blank No SIR (missing or invalid Final ID; MIC = ?, C, or X; unclaimed organism for Final ID)

* N indicates that the antimicrobial/organism does not have an upper breakpoint. When there is no upper breakpoint there are no criteria for calling an organism intermediate or resistant. This often occurs when there is an absence or rare occurrence of resistant strains of an organism. In this case, if the MIC is below the lower breakpoint the SIR results can be reported as susceptible but if the MIC is above the lower breakpoint the only result that can be reported is N or not susceptible.

I(nstrument SIR) does not appear if the EpiCenter system is attached and communicating. If communications with the EpiCenter system is lost, the field automatically appears.

(BD)E(Xpert SIR)

Possible field values are: Blank, S, I, R, N, and X. These values have the same meanings as I(nstrument SIR) above.

The BDXpert SIR field is not shown when the BDXpert System is disabled. It is also not shown for QC, ID, or ID/AST panels with AST disabled. This field contains the results calculated by the BDXpert System based on the execution of all enabled BDXpert Rules. Values appear in this field only if a BDXpert rule triggered and caused the BDXpert SIR value to differ from the instrument interpretation.

F(inal SIR)

Field values are the same as I(nstrument SIR) above. The Final SIR can be one of the following (from highest to lowest priority): a user-entered SIR value; the BDxpert SIR value; the instrument SIR value; or blank.

Final SIR values are produced when a MIC value or error is determined, AND there are no pending manual BDxpert rules (if the BDxpert System is enabled).

When the EpiCenter System is not connected and the BDxpert System is disabled, the Final SIR field is not displayed for an antimicrobial until the MIC value is something other than "?."

The chart below provides explanations of different combinations of blank and "X" results. Detailed explanations for actual results are provided as Special Messages (see Section 5.5.3).

Antimicrobial	MIC	Instrument SIR	Expert SIR		Example for Possible Causes
Any Drug	?	[Blank]	[Blank]	[Blank]	Results pending.
Any Drug	Phoenix MIC Result	S, I, R	[Blank]	S, I, R	Instrument SIR = Final SIR (No BDxpert rule SIR).
Any Drug	Phoenix MIC Result	S or I	R	R	Final SIR = BDxpert SIR if rule accepted.
Any Drug	Phoenix MIC Result	[Blank]	[Blank]	[Blank]	No breakpoints for this drug/organism combination within the chosen standard (CLSI, SFM, EUCAST, Custom).
Any Drug	Phoenix MIC Result	X	[Blank]	X	MIC value is outside the breakpoints for the selected standard. Example: Panel drug range = 1–16 µg/mL, susceptible breakpoint = 0.5 µg/mL
Any Drug	Phoenix MIC Result	X	[Blank]	X	SIR is suppressed by a BDxpert rule. User must provide Final SIR based on manual interpretation or additional testing.
Any Drug	X	[Blank]	[Blank]	[Blank]	The MIC for this drug and organism combination is not reported by the Phoenix system. An alternative method should be used.
All Drugs	X	[Blank]	[Blank]	[Blank]	This species is not included in the Phoenix AST taxonomy; perform an alternative method. An excessive amount of indicator was detected in the panel. The AST portion of the panel has been terminated and the isolate should be retested.
Any Drug	X	[Blank]	[Blank] or R	[Blank] or R	The MIC for this antibiotic is not reported (see Special Message). Interpretation based on BDxpert Rule.
Any Drug	C	[Blank]	R	R	Rapid Completion SIR (BDxpert SIR based on ID and/or another completed drug and/or resistance marker result).

Table 5-1 – MIC/SIR Values and Causes

Instrument Organism ID

Read-only field showing the organism identified by the instrument based on biochemical results in an ID or ID/AST panel. Organism IDs listed for a Tie (two possible organisms meet the results) are listed in the Instrument Organism ID field as "organism 1 / organism 2."

Confidence

The Confidence value computed by the system is based on the actual biochemical results versus the expected results. The Confidence value is a percentage from 0 to 99.

Supplemental Tests

Supplemental Tests are displayed in the Instrument Organism ID field and/or in the Special Messages display if there is more than one Organism listed there. Once you have performed these tests, the results of the tests will help to distinguish which organism ID to associate to the AST results. At this point, you may then select a single organism ID from the Final ID field.

Biochemical

Read-only field showing the abbreviation for the biochemical used to determine the ID.

Actual

Read-only field showing the observed biochemical result at the time when the organism ID was determined: + for positive or – for negative. A ? indicates that the biochemical result is pending or that the test was aborted prior to ID results being determined. X indicates an Error.

Expected

Read-only field showing the expected biochemical result according to the Instrument Organism ID: + or – for the organism. A V indicates that the result can be variable. This field is blank when the actual results are "?" (for non-QC panels), or until an Organism ID has been determined, or when more than one organism is listed in the Instrument Organism ID field.

5.5.1 BDXpert Triggered Rules Display

The BDXpert Triggered Rules display enables you to view the BDXpert system rules that have been triggered for a panel (the panel currently selected in Panel Results). You can view a listing of rules that have been triggered, the text of those rules, see the effect that the rules have on Final SIR values, and accept or reject pending (manual) rules. You can also re-execute all rules. See Figure 5–5.

After all rules have been reviewed and Accepted/Rejected, any changes must be saved in the Panel Results display via the "save" soft key.

The BDXpert Triggered Rules icon does not appear, and the display is not available if the EpiCenter system is attached and communicating with the Phoenix instrument. However, if communications with the EpiCenter system is lost, and the BDXpert system is reactivated in Instrument Configuration (Section 2.3.1), the icon reappears and the display can be accessed. When communications with the EpiCenter system is restored, BDXpert rules interpretations are once again performed at the EpiCenter system, and the "go to BDXpert Rules" icon no longer appears on the Panel Results display.

To access BDXpert Triggered Rules:

Press the “go to BDXpert Rules” soft key from the Panel Results screen

The display appears automatically if a change is made to the Final ID field (in Panel Results) that causes a manual BDXpert Rule to be triggered

BDXpert Triggered Rules soft keys:

BDXpert Triggered Rules fields:

Accession #

Read-only field that shows the current panel’s accession number.

Sequence #

Read-only field that shows the current panel’s barcode sequence number.

Isolate #

Read-only field that shows the current panel’s isolate number.

Final ID

Read-only field shows the panel’s final organism identification.

Rule Description Window

Read-only field that shows the text of the current rule. (The current rule is highlighted in the lower left window (Rules Window) of the display.) The color of the text of the rule indicates the following statuses:

GREEN	Automatic Rule
BLUE	Pending Manual Rule; execution of additional rules is suspended until you either Accept or Reject this rule
GREEN	Manual Rule that has been Accepted
RED	Manual Rule that has been Rejected

Rules Window

Rule

Read-only field that shows the numeric designation of the highlighted rule.
A maximum of 25 rules can be displayed in this window.

Status

Shows the status of the highlighted rule. The initial status of Automatic (a rule that executes automatically without user intervention) or Manual (a rule that must be manually Accepted or Rejected by the user) is set in the BDXpert Rule Configuration display (Section 2.3.3). Statuses shown in this window are:

Automatic – rule is enabled and set to Automatic

Pending – rule is enabled and set to Manual; Manual rules must be Accepted or Rejected; only the first Manual rule shows as Pending

Accepted – rule is enabled and set to Manual and has been Accepted by the user

Rejected – rule is enabled and set to Manual and has been Rejected by the user

Pending rules can be accepted or rejected via the "accept" or "cancel" key. Once a rule is accepted or rejected, the status can only be changed by re-executing the rules.

BDXpert TRIGGERED RULES			
Accession #	P010001	Isolate #	
Sequence #	411200000061	Final ID	B. pertussis
Rule Description Klebsiella oxytoca is usually susceptible to aminoglycosides. Verification of organism purity, identification and/or susceptibility is suggested.			
Rule #	Status	Antimier MIC	I F
1503	Automatic	GM >16	R
1522	Automatic	NM 16	R
		IPN 4	S
		CF >32	R
		FOX >16	R
		CTX >32	R
		FEP 16	I R

OK

Figure 5-5 – BDXpert Rules Triggered Display

Display reflects rules that had been triggered up to the time when the display was accessed. Display is not updated dynamically with new triggered rules.

AST Results Window

Antimicr(obial)

Read-only field that shows the abbreviation for the antimicrobial.

MIC

Read-only field that shows the MIC for that panel.

I(nstrument SIR)

Read-only field that shows the Instrument interpretation for that drug/MIC.

F(inal SIR)

Read-only field that shows the Final interpretation after any rules have been applied to the results.

5.5.2 Resistance Marker Display

The Resistance Marker display enables you to view any resistance markers that have triggered for a given panel. It is accessed through a soft key on the Panel Results display. See Figure 5 –6.

The Resistance Marker review display shows the same information in the top window as ID and AST Panel Results displays. In the bottom window, the following information is shown:

- The BDxpert rule number that triggered the Resistance Marker
- The Resistance Marker code (abbreviation)
- The Resistance Marker Name
- The BDxpert rule Description (green if the rule is automatic or accepted manual; blue if the rule is pending manual)

The resistance markers are shown in the order in which they trigger.

To access the Resistance Marker display:

Press the “toggle display” soft key from the Panel Results screen

The soft keys on the Resistance Marker display are the same as the Panel Results display.

0.5

PANEL RESULTS
CT01P

Accession #

Sequence #

Isolate #

Status

Final ID

Test Start

Test End

Location

Finalized ☐

Critical ☒

Rule #	RM	Name
1507	VRE	Vancomycin Resistant Enterococcus
0224	HLGR	High Level Gentamicin Resistant
0227	HLSR	High Level Streptomycin Resistant

Rule Description

For isolates with an interpretation of resistant for high levels of streptomycin (HLSR), synergy with cell wall active agents (e.g., ampicillin, penicillin, vancomycin) is not achievable.

Figure 5-6 – Resistance Marker Review Display

5.5.3 Special Messages Display

The Special Messages display enables you to access information about certain panel ID and AST results, as well as some panel readings. These special messages are triggered and are available for viewing regardless of whether the BDxpert System is enabled or disabled. If a recalled panel has an associated Special Message, the "special messages" soft key is displayed on the Panel Results display. See Figure 5-7.

Special Messages on a recalled panel are shown according to hierarchy on the Special Messages display, up to a maximum of 1000 characters. (Only complete special messages are shown.) If a message exceeds the 1000 character limit, it does not appear on the display. To view all messages in their entirety, print a Lab Report.

The Special Messages screen reflects messages that exist *at the time* the display is accessed: it is not updated dynamically with messages that are triggered after the display has been accessed. To view newly triggered messages, you must return to the Panel Results display, recall the panel again, and re-enter the Special Messages display again.

To access the Special Messages display:

- 1 From the Panel Results display, press the "more soft keys" soft key to view additional soft key definitions.

- 2 When the new soft keys appear, press the "special messages" soft key.

- 3 To return to Panel Results, press the "exit" soft key.


Figure 5-7 – Special Messages Display

5.5.4 Lab Report

The Lab Report contains all information for a panel sequence number that exists in the Phoenix database, including all information in the Panel Results screen, any special messages, BDxpert Rules that triggered, or Needs Attention Reasons if they exist. The Lab Report is not accessible from the Reports menu. It can only be printed from the Panel Results or Finalization displays.

The QC Lab Report is also accessible from the Panel Results display. It provides similar information to the Lab Report, but prints when a QC panel is being displayed and a report is requested.

The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), and Date and Time Printed, Software version (X.XXY) / Phoenix Update Disk (PUD) version (A.AAB)

Body of Report: Top Region: Accession #, Isolate #, Sequence #, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Lab Report is sorted by Accession # then by Isolate # within accessions.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDxpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDxpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report.

ID and AST sections of combination panels print on separate pages of the report.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

The system can be configured (General Instrument Configuration, Section 2.3.1.1) to print an "abbreviated" lab report. The abbreviated report does not contain the individual biochemical results for ID or ID/AST panels.

See Figure 5–8.

To print a Lab Report:

From the Panel Results display, press the "print" soft key.

Lab Report

03/30/03 09:11

VX.XXY / VA.AAB

Accession #

1313

Isolate #

1

Sequence #

420011234567

Panel Type

GPIDAST

Status

COMPLETE

Test Start

03/28/03 18:16

Test End

03/29/03 22:07

Instru # / Station

1/A01

Finalized

No

Final ID

Staphylococcus haemolyticus

Inoculum Density

0.5

Instrument ID(s)

Staph. haemolyticus

Confidence Value

99%

Biochemical

Instr

Expected

Result

Result

Biochemical

Instr

Expected

Result

Result

Biochemical

Instr

Expected

Result

Result

A_ARARR

•

•

A_GLPRB

•

•

A_LALT

-

•

A_LARGH

•

•

A_LHIST

•

•

A_LISO

•

•

A_LLEUH

•

•

A_LPHET

•

•

A_LPROB

•

•

•

•

•

•

•

[second page, header identical]

Antimicrobial

Instr

MIC

Instr

SIR

BDXpert

SIR

Final

SIR

Rule #

Ampicillin

>32

R

R

Azithromycin

>8

R

R

Ceftazidime

>64

R

R

Clindamycin

0.25

S

S

•

•

•

•

RESISTANCE MARKERS

abcd - ABRV Text of triggered resistance markers appears here.

BDXpert TRIGGERED RULES

1234 - AUTO Text of BDXpert triggered rule appears here.

End of Report

Figure 5–8 – Sample Lab Report

L-003342(R)

5-27

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

5.6 Finalization

The Finalization display enables you to review and finalize panel test results in a batch mode or individually. You can finalize panel results for one or all panels meeting the specified conditions, skip to a new panel, or modify and finalize panel results for a single panel. All fields on the Finalization display are read-only. To modify information for a panel, press the "panel results" soft key. When the instrument is connected to the EpiCenter Data Management System, the Finalization display is not available.

When you select Finalization, the instrument finds the records that are eligible for finalization. To be eligible, a panel must have a status of Removable, not be a QC panel, and have no unignored Needs Attention flags. Eligible panels are sorted first by accession number and then by isolate number.

When the "batch finalize" soft key is pressed, the hourglass icon indicates that the finalization operation is in progress. A maximum of 200 panels may be finalized in a given finalization session.

Two types of reports, standard Lab Report(s) and a finalization Summary Report, may be printed. The Lab Report can only be printed after one or more panels have been finalized. The Summary report may be printed any time.

See Figure 5-9.

To access Finalization:

Press the "finalization" soft key from the Main Status screen

Finalization soft keys:

"okay"

OK

"cancel"

"finalize all"

"finalize single"

"print Lab Report"

"print summary report"

"toggle display"

"panel results"

"more antimicrobials"

"exit"

Finalization fields:**Accession #**

Read-only field showing the panel's accession number.

Sequence #

Read-only field showing the panel's sequence number.

Isolate #

Read-only field showing the panel's isolate number.

Media

Read-only field showing the media type selected during Panel Login for Yeast ID panels.

Status

Read-only field showing the panel's status.

Final ID

Read-only field showing the panel's final organism ID.

Test Start

Read-only field showing the panel's start of testing date and time.

Test End

Read-only field showing the panel's end of testing date and time.

Location

Read-only field showing the panel's location.

Finalized

Read-only field showing the panel's finalized status.

Instrument Organism ID

Read-only field showing the instrument identification of the organism.

Confidence

Read-only field showing the confidence value for the Instrument Organism ID.

Biochemical

Read-only field showing the biochemical abbreviation.

Actual

Read-only field showing the actual biochemical result.

Expected

Read-only field showing the expected biochemical result.

Antimicrob(ial)

Read-only field showing the drug name.

MIC

Read-only field showing the Minimum Inhibitory Concentration value.


Figure 5–9 – Finalization Display

Interpretation

Read-only field showing the Susceptibility interpretation.

BDXpert SIR

Read-only field showing the final results calculated by the BDXpert System based on the execution of all enabled BDXpert Rules as if they had been executed automatically.

Final SIR

Read-only field showing the Final Susceptibility interpretation. To modify a Final SIR value, press the “panel results” soft key to go to the Panel Results display.

5.7 Configuration/Maintenance

Refer to Section 2.3 for complete information on system Configuration. Refer to Section 4.2 and Section 6 for information on system maintenance.

5.8 Reports

Numerous instrument reports are available for printing. The reports are as follows:

- ◆ Completed Lab Report – contains information for all panels whose status became “Complete” during the time period (up to the past 48 hours) you select.
- ◆ Accession Lab Report – provides information for a specified accession number.
- ◆ Needs Attention List Report – lists all panels in the instrument’s database that have a Needs Attention flag set that have not been ignored.
- ◆ Resident Panel Report – lists the panels contained in stations 1 – 25 for each tier.
- ◆ QC Lab Report – lists all Test Strain Organisms that have completed testing and all biochemical and/or antimicrobial MIC results (for a specified Panel Lot #) that exist in the Phoenix databases when the report is requested.
- ◆ Cumulative QC Report – provides information on quality control testing of all panel types.
- ◆ Daily Instrument Report – lists the status of the instrument at the time the report is generated, along with a maintenance checklist.
- ◆ Interpretation Rule Set Report – lists the antimicrobial breakpoints of the Interpretation Rule Set defined as the default Rule Set in the Instrument Configuration screen.
- ◆ BDXpert Rule Set Database Report – lists each BDXpert rule number and the text describing the rule, whether each rule is enabled/disabled and whether each rule shall trigger automatically/manually in the system.
- ◆ Organism ID Code List Report – displays all Organism Names, BD and LIS codes for all Organism Names that exist in the Phoenix database.
- ◆ Antimicrobial Code Report – displays all antimicrobials, BD and LIS codes for all antimicrobials, that exist in the Phoenix database from all panel configurations.
- ◆ Lab Report – contains all information for a panel sequence number that exists in the Phoenix database, including all information in the Panel Results screen, any special messages, BDXpert Rules that triggered, or Needs Attention Reasons if they exist. This report can only be printed from the Panel Results display or Finalization display.
- ◆ Finalization Summary Report – contains all the panels eligible for finalization at the time the report was requested, as well as finalization status. This report can only be printed from the Finalization display.
- ◆ Panel Inventory Lab Report – prints a Lab Report for all panels in the instrument. This report can only be printed from the Panel Inventory display.
- ◆ Custom Breakpoint Difference Report – prints a listing of differences between old breakpoints and new ones after a Phoenix Update Disk or Install/Upgrade operation. This report can only be printed from the Custom Interpretation Rule Set (Configuration) display.
- ◆ Current QC Panel Lot Report – contains information on the most recent QC test for each of the required strains for a panel lot, up to a maximum of 20 strains. This report can only be printed from the Panel Lot Definition display.
- ◆ Historical QC Panel Lot Report – contains information on all tests for a strain for the current instrument, up to 200 tests. This report can only be printed from the Panel Lot Definition display.
- ◆ Panel Lot Report – lists all the panel records for any panel lot number in the current instrument.
- ◆ Panel Lot Database Report – lists all the defined panel lots in the current instrument, and provides statistical and reference information on those lots.

Each of the reports is discussed in greater detail in the sections that follow.

5.8.1 Completed Lab Report

The Completed Lab Report contains information for *all panels* whose status became "Complete" during the time period (up to the past 48 hours) you select. The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), and Date and Time Printed, Software version (X.XXY) / Phoenix Update Disk (PUD) version (A.AAB).

Body of Report: Top Region: Accession #, Isolate #, Sequence #, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Completed Lab Report is sorted by Accession # then by Isolate # within accessions.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any Instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDxpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDxpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report. Any SIR values and Rule # are not reported for QC panels.

ID and AST sections of combination panels print on separate pages of the report.

See Figure 5–10.

To print a Completed Lab Report:

From the Main Status screen, press the "print" soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Completed Lab Report.

Start Date/Time

This field defaults to the current date and time. This is the *end time* of the Completed Lab Report print period. The beginning time is the hour offset specified in Prior Period (below). Press the ▲ or ▼ key to increase or decrease the displayed value (Month, Day, Year, Hour, or Minute) to the desired end date/time.

Prior Period

This field represents the number of hours back to print the report. That is, if you select "2" for this field, the system prints Lab Reports for all panels that have completed testing from the time specified in **Start Date/Time back 2 hours**. This field defaults to 24 hours. Press the ▲ or ▼ key to increase or decrease the displayed value until the desired number of hours is shown. You can select from 1 to 48 hours ago.

QC Checkbox

Press "select" to print a Completed QC Lab Report (disabled by default).

Press the "print" soft key again to print the report.

Completed Lab Report

03/30/03 09:11
VX.XXY / VA.AAB

Accession #1313

Isolate #1

Sequence #420011234567

Panel TypeGPIDAST

StatusCOMPLETE

Test Start03/28/03 18:16

Test End03/29/03 22:07

Instru # / Station1/A01

FinalizedNo

Final IDStaphylococcus haemolyticus

Inoculum Density0.5

Instrument ID(s)Staph. haemolyticus

Confidence Value99%

Supplemental Test(s)

Biochemical

Instr

Expected

Result

Result

A_ARARR

-

-

A_GLPRB

-

-

Biochemical

Instr

Expected

Result

Result

A_LARGH

-

-

A_LHIST

-

-

Biochemical

Instr

Expected

Result

Result

A_LLEUH

-

-

A_LPHET

-

-

A_LALT

-

-

A_LISO

-

-

A_LPROB

-

-

.

.

.

[second page, header identical]

Antimicrobial

Instr

MIC

Ampicillin

>32

Instr

SIR

Azithromycin

R

BDXpert

SIR

Ceftazidime

R

Final

SIR

Clindamycin

R

Rule #

.

RESISTANCE MARKERS

abcd - ABRV Text of triggered resistance markers appears here.

BDXpert TRIGGERED RULES

1234 - AUTO Text of BDXpert triggered rule appears here.

End of Report

Figure 5-10 – Sample Completed Lab Report

5-34

L-003342(R)

ООО ЦРПИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

5.8.2 Accession Lab Report

The Accession Lab Report is basically a collection of lab reports for a specified accession number. It provides information for a specified accession number. Such information includes, as applicable: organism ID results, including specific biochemical reactions; AST results including SIR interpretation and MIC; QC pass/fail results; and any BDxpert Rules that were triggered. The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), Date and Time Printed, Software version (X.XXY) / Phoenix Update Disk (PUD) version (A.AAB).

Body of Report: Top Region: Accession #, Isolate #, Sequence #, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Completed Lab Report is sorted by Accession # then by Isolate # within accessions.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any Instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDxpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDxpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report. Any SIR values and Rule # are not reported for QC panels.

ID and AST sections of combination panels print on separate pages of the report.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

See Figure 5-11.

To print an Accession Lab Report:

From the Main Status screen, press the "print" soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Accession Lab Report.

Accession #

Type in or scan the desired Accession number. If you enter an incorrect number, you can press the "clear" soft key to clear the field and reenter the correct number.

Press the "print" soft key again to print the report.

Accession Lab Report

03/30/03 09:11
VX.XXY / VA.AAB

Accession #
Isolate #
Sequence #
Panel Type
Status

1313
1
420011234567
GPIDAST
COMPLETE

Test Start
Test End
Instru # / Station
Finalized

03/28/03 18:16
03/29/03 22:07
1/A01
No

Final ID

Staphylococcus haemolyticus

Inoculum Density

0.5

Instrument ID(s)
Staph. haemolyticus

Confidence Value
99%

Biochemical	Instr Result	Expected Result	Biochemical	Instr Result	Expected Result	Biochemical	Instr Result	Expected Result
A_ARARR	•	•	A_GLPRB	•	•	A_LALT	•	•
A_LARGH	•	•	A_LHIST	•	•	A_LISO	•	•
A_LLEUH	•	•	A_LPHET	•	•	A_LPROB	•	•
.			.			.		
.			.			.		

[second page, header identical]

Antimicrobial	Instr MIC	Instr SIR	BDXpert SIR	Final SIR	Rule #
Ampicillin	>32	R		R	
Azithromycin	>8	R		R	
Ceftazidime	>64	R		R	
Clindamycin	0.25	S		S	
.		.			
.		.			

RESISTANCE MARKERS

abcd - ABRV Text of triggered resistance markers appears here.

BDXpert TRIGGERED RULES

1234 - AUTO Text of BDXpert triggered rule appears here.

End of Report

Figure 5-11 – Sample Accession Lab Report

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5.8.3 Needs Attention List Report

The Needs Attention List Report lists all the panels in the instrument's database that have a Needs Attention flag set that have not been ignored. This report can also be printed from the Needs Attention display. The report provides the following information:

Header: Report Title, Filtered report notification (if selected), Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Needs Attention Reason, Sequence #, Accession #, Isolate #, Instr #/Station (location), and Status (Ongoing, Complete). If the report is filtered (default selection), an asterisk appears to the left of the Reason for panels with multiple Needs Attention conditions. The report is sorted by the priority of the Needs Attention Reasons (see Section 5.11), and by Accession within each Reason type.

See Figure 5-12.

To print a Needs Attention List Report:

From the Main Status screen, press the "print" soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Needs Attention List Report.

If you want to print an unfiltered report, press **Tab** to advance to the Filtered checkbox. Press **Spacebar** or the "select" soft key to clear the check in the box.

Press the "print" soft key again to print the report.

Needs Attention List Report					
Filtered on highest priority reason				10/27/03 09:11	
Reason	Sequence #	Accession #	Isolate #	Instr # / Station	Status
No Growth On Panel	420011234567	1313	1	1/A02	COMPLETE
Review QC Results	420027654321	1344	1	1/B03	COMPLETE
Can Not Determine Organism ID	420031234123	1489	1	1/C04	ONGOING
Can Not Identify Barcode	420047654765	1566	1	1/D05	COMPLETE
End of Report					

Figure 5–12 – Sample Needs Attention List Report

5.8.4 Resident Panel Report

The Resident Panel Report lists the panels contained in stations 1 – 25 for each tier detected during the last inventory scan. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument #.

Body of Report: Accession #, Isolate #, Sequence #, QC (if panel is QC), Test Start with time, Inoculum Density, Panel Type, Status (Ongoing, Complete), and the highest priority Needs Attention Reason if one exists. The report is sorted by Accession # and then by Isolate # within each Accession.

See Figure 5–13.

To print a Resident Panel Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Resident Panel Report.

Press the “print” soft key again to print the report.

Resident Panel Report								
07/30/99 09:11								
Instrument #1								
Accession #	Isolate #	Sequence #	QC	Test Start	Inoculum Density	Panel Type	Status	Needs Attention
1654	1	420011234567		07/28/99 14:52	0.5	GNIDAST	COMPLETE	Missing Accession Number
1876	1	420011234568		07/28/99 15:33	0.5	GNIDAST	COMPLETE	.
11899	1	420011234569		07/28/99 15:36	0.5	GNIDAST	COMPLETE	.
1945	1	420011234560		07/28/99 15:45	0.5	GNIDAST	COMPLETE	.
.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.
End of Report								

Figure 5-13 – Sample Resident Panel Report

5.8.5 QC Lab Report

The QC Lab Report lists all QC panels from the Test Start date field entry to the current date. It lists all Test Strain Organisms that have completed testing and all biochemical and/or antimicrobial MIC results (for a specified panel lot number) that exist in the Phoenix database. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed, Software version (X.XXY) / Phoenix Update Disk (PUD) version (A.AAB).

Body of Report: Panel Lot # and Expiration Date, Test Start and time, Test End and time, Panel Type, Instr #/Station (location), Status (Ongoing, Complete), Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, Sequence #, Accession #, Isolate #, Test Strain, Inoculum Density, Media Type (Yeast ID panels only), Instrument ID(s), and QC Status of PASS/FAIL/REPEAT/REVIEW. At the bottom of the report, any Needs Attention reasons or Special Messages are printed. Each Biochemical, along with Instr(ument) Result, and Expected Result are provided, as well as Antimicrobials, Instr(ument) MICs, Expected MICs, and Pass/Fail status.

This report is only available when BD EpiCenter is disabled.

See Figure 5–14.

To print a QC Lab Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight QC Lab Report.

Panel Lot #

Scan or type in the Panel Lot #; if QC Lot Support is enabled, press TAB to advance to the Panel Lot # field, then press DOWN ARROW to drop down the list of all defined lot numbers, and highlight the desired lot number then press the “select” soft key to select the lot.

Test Start

In the first Test Start date field, press UP ARROW or DOWN ARROW to enter the desired day, month, or year (depending on your configuration). Press  to advance to the next date field and press UP ARROW or DOWN ARROW to enter the desired day, month, or year. Press  to advance to the next date field, and press UP ARROW or DOWN ARROW to enter the final portion of the date. If you enter an incorrect number, you can press the “clear” soft key to clear all the fields and reenter the correct number.

If you do not enter a Test Start date, the instrument prints all data from the database for the Panel Lot # entered.

Press the “print” soft key again to print the report.

QC Lab Report

10/30/03 09:11
VX.XXY / VA.AAB

Panel Lot #	777777	Expiration Date	10/07/05
Test Start	10/28/03 10:02	Test End	10/29/03 22:07
Panel Type	GPIDAST	Instru # / Station	1/A01
Status	COMPLETE	Tech ID	ljk
ID Broth Lot #	6666666	Expiration Date	07/07/05
AST Broth Lot #	5555555	Expiration Date	08/08/05
Indicator Lot #	4444444	Expiration Date	09/09/05
Sequence #	42001222222		
Accession #	9911		
Isolate #	1		

Test Strain 29213 Staphylococcus aureus **Inoculum Density** 0.5

Instrument ID(s)
Staph. aureus

QC Status PASS

<u>Biochemical</u>	<u>Instr</u>	<u>Expected</u>	<u>Biochemical</u>	<u>Instr</u>	<u>Expected</u>	<u>Biochemical</u>	<u>Instr</u>	<u>Expected</u>
	<u>Result</u>	<u>Result</u>		<u>Result</u>	<u>Result</u>		<u>Result</u>	<u>Result</u>
A_ARARR	-	•	A_GLPRB	•	•	A_LALT	•	•
A_LARGH	•	v	A_LHIST	•	•	A_LISO	•	•
A_LLEUH	•	•	A_LPHET	•	-	A_LPROB	-	•

[second page, header identical]

<u>Antimicrobial</u>	<u>Instr</u>	<u>Expected</u>	<u>Pass/Fail</u>
	<u>MIC</u>	<u>MIC</u>	
Amikacin	2	1.00-4.00	Pass
Chloramphenicol	4	2.00-8.00	Pass
Ciprofloxacin	0.25	0.13-0.50	Pass
Cephalexin	4	2.00-8.00	Pass

Needs Attention Reasons

Text of Needs Attention Reasons goes here

SPECIAL MESSAGES

The beta-lactamase test (nitrocefin based) passed QC.

End of Report

Figure 5-14 – Sample QC Lab Report

5.8.6 Cumulative QC Report

The Cumulative QC Report provides information on completed quality control testing of all panel types. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument #.

Body of Report: Selection Criteria: Panel Lot #, Panel Type, and Test Strains selected.

Below this, the panel Sequence # (sort order), QC Status (PASS, FAIL, REVIEW, ERROR), Test Strain, Test Start and time, Panel Lot #, ID Broth Lot #, AST Broth Lot #, Indicator Lot #, and Tech ID are listed for each panel.

This report is only available when BD EpiCenter is disabled.

See Figure 5–15.

To print a Cumulative QC Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Cumulative QC Report.

If QC Lot Support is disabled:

Panel Lot #	Type in or scan the desired panel lot number. If you enter an incorrect number, you can press the “clear” soft key to clear all the fields and reenter the correct number.
All Checkbox	Initially, this box is checked, indicating that all completed QC data will be printed. Press the “move to next field” soft key to advance the cursor to this checkbox and press the “select” soft key to UN-check the box to specify panel types or test strains on which to report. The All checkbox supercedes any entry in the Panel Lot # field.
Panel Type	Press UP ARROW or DOWN ARROW to highlight the desired panel type, or select “(all)” for all types. Only completed panel types appear in the list.
Test Strain	Press UP ARROW or DOWN ARROW to highlight the desired test strain, or select “(all)” for all strains.

If QC Lot Support is enabled:

Panel Lot #	To select results to print by lot number, use this field. Select “(all)” to print all panel lot numbers. To specify a specific lot #, press DOWN ARROW to drop down the list of defined lot numbers, highlight the desired lot number then press the “select” soft key to select the lot. The panel type field is completed automatically.
Panel Type	To select results to print by panel type, use this field. Press DOWN ARROW to select a specific panel type to print, or select “(all)” for all types.

Test Strain

To select results to print by test strain, use this field. Press DOWN ARROW to highlight the desired test strain, or select "(all)" for all strains. The "(all)" selection prints all strains when a panel lot number is specified (even strains not required for a panel type).

Press the "print" soft key again to print the report.

Cumulative QC Report								
10/30/03 09:11								
Instrument #1								
Panel Lot #	(all)							
Panel Type	(all)							
Test Strain	(all)							
Sequence #	QC Status	Test Strain	Test Start	Panel Lot #	ID Broth Lot #	AST Broth Lot #	Indicator Lot #	Tech ID
420011234567	PASS	ATCC51299	10/28/03 14:52	1111111	5555555	9999999	7878787	ljk
420011234568	PASS	ATCC25923	10/28/03 15:33	2222222	6666666	1212121	9090909	bec
420011234569	PASS	ATCC29213	10/28/03 15:36	3333333	7777777	3434343	1234567	me
420011234560	PASS	ATCC29062	10/28/03 15:45	4444444	8888888	5656565	2345678	up
.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.
End of Report								

Figure 5-15 – Sample Cumulative QC Report

5.8.7 Daily Instrument Report

The Daily Instrument Report lists the status of the instrument at the time the report is generated, and provides areas to record maintenance activities. The Daily Instrument Report can be set to print automatically at a time you specify. Refer to Section 2.3.1.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Software version (X.XXY) / Phoenix Update Disk (PUD) version (A.AAB).

Body of Report: Instrument #, Serial #, Instrument Temperature Pass/Fail status, Carousel Rotational Test Pass/Fail status, Power Supply Check Pass/Fail status, Normalizer Panels Sequence #, Pass/Fail status, and Expiration Status (date if expiration is more than 60 days; "expires on date" if expiration is between 60 and 0 days; and "EXPIRED" if the panel is expired) for each tier, and blanks to record the reading, Pass/Fail status, and Tech ID for each of the following maintenance checks: Daily: Instrument Temperature (Main Status screen), Standard Panel Temperature, Printer Paper Supply; Weekly: Internal Green LEDs, Internal Red LEDs, Internal Amber LEDs, System Alert Indicator, and Instrument Audible Alarm. An area is provided for comments at the bottom.

The instrument temperature is considered to have passed when there are no outstanding E01 temperature alerts.

See Figure 5–16.

To print a Daily Instrument Report:

From the Main Status screen, press the "print" soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Daily Instrument Report.

Press the "print" soft key again to print the report.

Daily Instrument Report

11/30/07 09:11
VX.XXY / VA.AAB

Instrument #1
Serial #12345

Instrument TemperaturePass
Carousel Rotational TestPass
Power Supply CheckPass

Normalizer Panels

Tier	Sequence #	Status	Expiration Status
A	429990000001	Pass	10/21/08
B	429990000002	Pass	10/21/08
C	429990000003	Pass	10/21/08
D	429990000004	Pass	10/21/08

Daily	Pass	Fail	Tech ID
Record Instrument Temperature From Main Screen Range: (33.5 - 36.5)			
Record Standard Panel Temperature Range: (33.5 - 36.5)			
Check Printer Paper Supply			
Weekly			
Instrument Internal Green LEDs			
Instrument Internal Red LEDs			
Instrument Internal Amber LEDs			
Instrument System Alert Indicator			
Instrument Audible Alarm			
Comments			

End of Report

Figure 5-16 – Sample Daily Instrument Report

5.8.8 Interpretation Rule Set Report

The Interpretation Rule Set Report lists the antimicrobial breakpoints of the currently selected Interpretation Rule Set (defined as the default Rule Set in the Instrument Configuration screen). The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Rule Set, Rule Version, and columns for Antimicrobial (sort order), Test Group, Organism Group, Organism Name, S(usceptible) value, and R(esistance) value. Each antimicrobial breakpoint is listed in a separate row of the report.

The Interpretation Rule Set is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

To print reports of ALL the rule sets in the instrument, first select a rule set in Instrument Configuration (Section 2.3.1), then access the Reports menu and select Interpretation Rule Set Report. The currently selected rule set prints. When printing is complete, return to Instrument Configuration, and select the next rule set. Then access the Reports menu again and print the current Interpretation Rule Set. Continue selecting rule sets and printing until all rule set selections are printed. Remember to return to Instrument Configuration and select the desired rule set to use for interpretations when all printing is complete.

NOTE

Modification of the Interpretation Rule Set should not be performed while there are ongoing panels. This could lead to inaccurate interpretations.

See Figure 5–17.

To print an Interpretation Rule Set Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Interpretation Rule Set Report.

Press the “print” soft key again.

OK

A message informs you about the report size. Press the “okay” soft key to print the report, or press the “cancel” soft key to cancel printing.

Interpretation Rule Set Report					
03/30/05 09:11					
Rule Set	CLSI				
Rule Version	M100-S13				
Antimicrobial	Test Group	Organism Group	Organism Name	S	R
Ampicillin	A	ENTC_IC	(all)	<=8.00	>=16.00
Ampicillin	A	ENTERIC_IC	(all)	<=8.00	>=32.00
Ampicillin	A	HAEINF_IC	(all)	<=1.00	>=4.00
Ampicillin	A	HAEIOTHF_IC	(all)	<=1.00	>=4.00
Ampicillin	A	LISMON_IC	(all)	<=2.00	-
Ampicillin	O	STAAUE_IC	(all)	<=0.25	>=0.50
Ampicillin	O	STAOTH_IC	(all)	<=0.25	>=0.50
Ampicillin	A	STRBET_IC	(all)	<=0.25	-
.	.	.	.	.	.
.	.	.	.	.	.
.	.	.	.	.	.
End of Report					

Figure 5-17 – Sample Interpretation Rule Set Report (sample data only)

Organism (Interpretation) Group Codes:
CLSI, SFM, or EUCAST

Interpretation Code	Interpretation Name
ACIN_IC	<i>Acinetobacter</i> spp.
AERM_IC	<i>Aeromonas</i> spp.
AERMHC_IC	<i>Aeromonas hydrophila</i> complex
BACI_IC	<i>Bacillus</i> spp.
BURCEP_IC	<i>Burkholderia cepacia</i>
CARHOM_IC	<i>Cardiobacterium hominis</i>
CORY_IC	<i>Corynebacterium</i> spp.
EIKCOR_IC	<i>Eikenella corrodens</i>
ENTC_IC	<i>Enterococcus</i> spp.
ENTERIC_IC	Enterobacteriaceae
ERYRHU_IC	<i>Erysipelothrix rhusiopathiae</i>
GNCBOTH_IC	Gram-negative coccobacilli, including <i>Moraxella</i> spp. other than MORBCAT_IC
GNROTH_IC	Gram-negative rod, other than ACIN_IC, BURCEP_IC, ENTERIC, NFGNROTH_IC, PSEAER_IC, STEMAL_IC, PAST_IC
GPCOTH_IC	Gram-positive, other than , <i>Staphylococcus</i> spp., <i>Streptococcus</i> spp., ENTC_IC, LEUC_IC, PEDI_IC
GPROTH_IC	Gram-positive rod, other than BACI_IC, CORY_IC, ERYRHU_IC
KIN_IC	<i>Kingella</i> spp.
LEUC_IC	<i>Leuconostoc</i> spp.
LISMON_IC	<i>Listeria monocytogenes</i>
LISTOTH_IC	<i>Listeria</i> spp., other than LISMON_IC
MORBCAT_IC	<i>Moraxella (Branhamella) catarrhalis</i>
NFGNROTH_IC	Nonfermentative GNR, other than ACIN_IC, BURCEP_IC, PSEAER_IC, STEMAL_IC, ACTBACT_IC, CARHOM_IC, EIKCOR_IC
NO_INTERP_IC	No Interpretation
PAST_IC	<i>Pasteurella</i> spp.
PEDI_IC	<i>Pediococcus</i> spp.
PLESI_IC	<i>Plesiomonas shigelloides</i>
PSEAER_IC	<i>Pseudomonas aeruginosa</i>
STAAUE_IC	<i>Staphylococcus aureus</i>
STAOTH_IC	<i>Staphylococcus</i> spp., other than STAAUE_IC
STEMAL_IC	<i>Stenotrophomonas maltophilia</i>
STRBET_IC	<i>Streptococcus</i> beta-hemolytic
STROTH_IC	<i>Streptococcus</i> spp., other than STRBET_IC, STRPNE_IC, STRVIR_IC
STRPNE_IC	<i>Streptococcus pneumoniae</i>
STRVIR_IC	<i>Streptococcus viridans</i> group
VIBCHO_IC	<i>Vibrio cholerae</i>
VIBROTH_IC	<i>Vibrio</i> spp., other than VIBCHO_IC

Table 5-2 – Organism Group Codes

5.8.9 BDXPert Rule Set Database Report

The BDXPert Rule Set Database Report lists each BDXPert rule number and the text describing the rule, whether each rule is enabled/disabled and whether each rule shall trigger automatically/manually in the system. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Rule Set (CLSI, SFM, EUCAST, or Custom), and Based On (CLSI, SFM, or EUCAST if Rule Set is Custom).

Body of Report: Rule #, text of the rule, Enabled/Disabled status, and Automatic/Manual status.

The BDXPert Rule Set Database is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

See Figure 5–18.

To print a BDXPert Rule Set Database Report:

- 

From the Main Status screen, press the “print” soft key to access the Reports menu.
- 

Use the UP ARROW or DOWN ARROW to highlight BDXPert Rule Set Database Report.
- 

Press the “print” soft key again.
- OK

A message informs you about the report size. Press the “okay” soft key to print the report, or press the “cancel” soft key to cancel printing.

BDXpert Rule Set Database Report			
Rule Set Based On		CLSI	
Rule #		03/17/05 09:11	
101	Only results for ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported for fecal isolates of Salmonella spp. or Shigella spp. No other drugs should be reported.	Enabled	Automatic
103	Results for chloramphenicol and a third generation cephalosporin should be reported for Salmonella spp. from extraintestinal sites.	Enabled	Automatic
104	Salmonella spp. and Shigella spp. should not be reported as susceptible or intermediate to first and second generation cephalosporins even if they appear to be susceptible in vitro because they are not effective clinically.	Enabled	Automatic
105	Ampicillin is the class representative for amoxicillin.	Enabled	Automatic
106	.		
.			
.			
.			
End of Report			

Figure 5–18 – Sample BDXpert Rule Set Database Report

5.8.10 Organism ID Code List Report

The Organism ID Code List Report prints all Organism Names and Abbreviations for all Organism Names that exist in the Phoenix database. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Organism name (sort order), BD Code (abbreviation), and LIS Code (if enabled).

The Organism ID Code List is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

See Figure 5–19.

To print an Organism ID Code List Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Organism ID Code List Report.

Press the “print” soft key again.

OK

A message informs you about the report size. Press the “okay” soft key to print the report, or press the “cancel” soft key to cancel printing.

Organism ID Code List Report		
07/30/99 09:11		
<u>Organism</u>	<u>BD Code</u>	<u>LIS Code</u>
Achr. species	ACHRSPE	
Acinet. baumannii	ACINBAU	
Acinet. baumannii/calco. cplx	ACINBCX	
Acinet. calcoaceticus	ACINCAL	
.	.	
.	.	
.	.	
End of Report		

Figure 5–19 – Sample Organism ID Code List Report

5.8.11 Antimicrobial Code Report

The Antimicrobial Code Report prints all antimicrobials and abbreviations for all antimicrobials that exist in the Phoenix database from all panel configurations. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Antimicrobial name (sort order), BD Code (abbreviation), and LIS Code (if enabled).

The Antimicrobial Code is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

See Figure 5–20.

To print an Antimicrobial Code Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Antimicrobial Code Report.

Press the “print” soft key again.

OK

A message informs you about the report size. Press the “okay” soft key to print the report, or press the “cancel” soft key to cancel printing.

Antimicrobial Code Report		
03/30/05 09:11		
Antimicrobial	BD Code	LIS Code
Amikacin	AN	AN
Amoxicillin	AMX	AMX
Amoxicillin/Clavulanate	AMC	AMC
.	.	.
.	.	.
.	.	.
End of Report		

Figure 5–20 – Sample Antimicrobial Code Report

5.8.12 Lab Report/QC Lab Report

The Lab Report (Section 5.5.4) contains all information for a panel sequence number that exists in the Phoenix database, including all information in the Panel Results screen, any special messages, BDxpert Rules that triggered, or Needs Attention Reasons if they exist. The Lab Report is not accessible from the Reports menu. It can only be printed from the Panel Results or Finalization displays. A Panel Inventory Lab Report prints the same information for all panels listed in the Panel Inventory display (i.e., resident in the instrument and with final or partial results). Refer to Section 5.5.4 for additional information.

The QC Lab Report (Section 5.8.5) is also accessible from the Panel Results display. It provides similar information to the Lab Report, but for QC panels.

5.8.13 Finalization Summary Report

The Finalization Summary Report contains a listing of all the panels eligible for finalization at the time the report was requested, as well as finalization status. The Finalization Summary Report is not accessible from the reports menu. It can only be printed from the Finalization display. Refer to Section 5.6 for additional information.

5.8.14 Custom Breakpoint Difference Report

The Custom Breakpoint Difference Report contains a listing of differences between old breakpoints and new ones after a Phoenix Update Disk or install/upgrade operation. The Custom Breakpoint Difference Report is not accessible from the Reports menu. It can only be printed from the Custom Interpretation Rule Set (Configuration) display. Refer to Section 2.3.3 for additional information.

5.8.15 Current QC Panel Lot Report

The Current QC Panel Lot Report contains information on the most recent QC test for each of the required strains for a panel lot, up to a maximum of 20 strains. The report includes information for any instruments whose data has been restored to the current instrument. The Current QC Panel Lot Report cannot be printed from the Reports menu. It can only be printed from the Panel Lot Definition display. Refer to **Section 2.3.6 – Panel Lot Definition**.

5.8.16 Historical QC Panel Lot Report

The Historical QC Panel Lot Report contains information on all tests for a strain for the current instrument (only), up to 200 tests. The Historical QC Panel Lot Report cannot be printed from the Reports menu. It can only be printed from the Panel Lot Definition display. Refer to **Section 2.3.6 – Panel Lot Definition**.

5.8.17 Panel Lot Report

The Panel Lot Report lists all the panel records for any panel lot number in the current instrument. The report first lists clinical panels, then QC panels. Within each of those groups, the report is sorted by Accession # then Isolate #.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Instrument where printed.

Body of Report: Panel Lot # and Panel Type; Accession #, Isolate #, Sequence #, QC (if panel is QC), Test Date, and Status (Pending, Ongoing, Complete) for each panel tested that belongs to the lot.

This report is only available when BD EpiCenter is disabled and QC Lot Support is enabled.

See Figure 5–21.

To print a Panel Lot Report:

From the Main Status screen, press the "print" soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Panel Lot Report.

Press the TAB key to move to the Panel Lot field.

Press DOWN ARROW to drop down the list of Panel lots; then highlight the desired lot and press the "select" soft key to select the highlighted lot.

Press the "print" soft key again.


Figure 5-21 – Sample Panel Lot Report

5.8.18 Panel Lot Database Report

The Panel Lot Database Report lists all the defined panel lots in the current instrument, and provides statistical and reference information on those lots.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Instrument where printed.

Body of Report: Panel Lot # (sort order, descending); Panel Type; Expiration Date; Extension Date (if Expiration date was extended); Start and End Sequence #s (Range); Definition Date; First and Last Date Used; and number of Panels Used.

This report is only available when BD EpiCenter is disabled and QC Lot Support is enabled.

See Figure 5–22.

To print a Panel Lot Database Report:

From the Main Status screen, press the “print” soft key to access the Reports menu.

Use the UP ARROW or DOWN ARROW to highlight Panel Lot Database Report.

Press the “print” soft key again.

Panel Lot Database Report										
07/15/04 09:11										
Instrument #1										
Panel Lot #	Panel Type	Expiration Date	Extension Date	Sequence Number Range		Definition Date	First Used Date	Last Used Date		Panels Used
				Start	End					
9856543	NMIC/ID-1	08/14/05	.	425010000001	425010000099	06/12/04	06/22/04	07/14/04	.	5
7600001	SMIC/ID-2	04/01/06	.	428510004200	428510004210	06/12/04	06/22/04	07/14/04	.	2
7500001	PMIC/ID-17	04/01/05	.	425350000101	425350000200	06/12/04	06/23/04	07/14/04	.	5
7400001	NMIC	10/06/08	.	422509999990	422509999999	06/12/04	06/23/04	07/14/04	.	3
.	.	.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.	.	.
.	.	.	.	.	.	.	.	.	.	.
End of Report										

Figure 5–22 – Sample Panel Lot Database Report

5.9 Panel Inventory

The Panel Inventory display enables you to view a listing of all panels in the instrument (except temperature reference panels). This listing can be sorted in ascending (0 – 9, A – Z) or descending (9 – 0, Z – A) order by the following fields: panel sequence number, accession/isolate number, results, or needs attention reason. The initial default sort order is by accession/isolate number in ascending order. Subsequently, the list defaults to the last sort criteria and display configuration (primary/secondary) you used. If there are no panels in the instrument, the display shows the message, “No Data Available.”

The top of the Panel Inventory display shows the number of panels in the list (if there is more than one panel), the sort field, and whether the sort field is ascending or descending. See Figure 5–23.

The following fields are shown on the Panel Inventory display (if the information is known):

- ①② Sequence # of the panel
- ①② Accession # (panels without accession numbers are listed first)
- ①② Iso(late number)
- ① (Inoculum) Density (blank for AST panels; “?” for ID panels until first test completes)
- ① Status (ongoing; complete; rapid if Rapid Completion is enabled and has been triggered for a panel)
- ①② Results (**final** if panel is complete and there are no active Needs Attention conditions [all MIC results are determined for an AST panel, or all MIC values *and* the organism ID are determined for a Combination panel, or the organism ID is determined for an ID panel]; **partial** if a panel is ongoing or complete but has an active (unignored) Needs Attention condition [at least one MIC value is determined for an AST or Combination panel, or the organism ID is determined for a Combination panel]; **none** if no MIC values or organism ID is determined for any type of panel; a + preceding the result indicates that a Resistance Marker has triggered for that panel)
- ① N(eeds) A(ttention)/C(ritical Panel) (“*” if a Needs Attention reason that has not been ignored exists; **blank** when no Needs Attention reasons exist or they have all been ignored; ♦ if a critical panel event is triggered)
- ② Final (Organism) ID

Fields listed above with a “one” (①) appear on the primary Panel Inventory display. Fields listed above with a “two” (②) appear on the secondary display (Figure 5–24), which is accessed from the primary display via a soft key.

After the list is displayed, you can use the ▲ or ▼ key to highlight a panel and access the Panel Results display to view or edit panel information, to perform an instrument “locate panel” operation, or to print a Lab Report for panels with final or partial results. These functions are performed by the soft keys described below. Panels shown in red are critical panels that have not been acknowledged.

PANEL INVENTORY						
15 Panels Sorted By Accession # Ascending						
Sequence #	Accession #	Iso	Density	Status	Results	NA/C
425810003591	3952		0.5	Complete	+Final	
425810003588	A49619	4	0.5	Complete	+Final	
425190008038	EN11002	6	0.5	Complete	+Final	
425190008902	EN11021	6	0.5	Complete	+Partial	*
425200001443	QCA29213	4	0.5	Complete	+Final	
425200001452	QCA43300	6	0.5	Complete	+Partial	*
425200001468	QCA51299	6	0.5	Complete	+Final	
428510004208	41583	1	0.5	Complete	Final	
425190008891	EN10909	4	0.5	Complete	Partial	*
425190000025	EN10909	4	0.5	Complete	Partial	*
425190008898	EN11764	7	0.5	Complete	Partial	*
425190008900	EN11807	7	0.5	Complete	Partial	*
425190008021	QCA25922	6	0.5	Complete	Final	
425190008910	SEC2009	4	0.5	Complete	Partial	*
425190008899	SEFR197	6	0.5	Complete	Final	

Figure 5-23 – Panel Inventory Display (primary)

PANEL INVENTORY				
15 Panels Sorted By Accession # Ascending				
Sequence #	Accession #	Iso	Final ID	Results
425810003591	3952	4	Streptoc. pneumoniae	+Final
425810003588	A49619	4	Streptoc. pneumoniae	+Final
425190008038	EN11002	6	Esch. coli	+Final
425190008902	EN11021	6	Kleb. pneumo. ssp pneu.	+Partial
425200001443	QCA29213	4	Staph. aureus	+Final
425200001452	QCA43300	6	Staph. aureus	+Partial
425200001468	QCA51299	6	Enteroc. faecalis	+Final
428510004208	41583	1	Streptoc. mitis	Final
425190008891	EN10909	4	Kleb. oxytoca	Partial
425190000025	EN10909	4	Kleb. oxytoca	Partial
425190008898	EN11764	7	Salm. species	Partial
425190008900	EN11807	7	Cit. farmeri	Partial
425190008021	QCA25922	6	Esch. coli	Final
425190008910	SEC2009	4	Pseud. aeruginosa	Partial
425190008899	SEFR197	6	Ser. marcescens	Final

Figure 5-24 – Panel Inventory Display (secondary)

To access Panel Inventory:

Press the "panel inventory" soft key from the Main Status screen

Panel Inventory soft keys:

"re-execute"

"results"

"locate panel"


"print Lab Report"

"toggle display"

"sort ascending/
descending"

"toggle keys"


"sort by sequence #"

"sort by accession #"


"sort by results"

"sort by needs attention"

"exit"

5.10 System Alerts

The System Alerts display enables you to review any existing system alerts that may have occurred or that may still exist in the instrument.

When you access the System Alerts display, the first alert in the list is shown in the middle of the screen. This alert is represented by an "E" code (such as E12). Below the "Enn" (where nn is a two-digit number) code one or more subcodes appear which are 8 digits. The subcode represents the specific condition of the alert. For example, E07 is a power supply alert, and subcode 00000002 specifies a problem in the 5 Volt power supply.

To view additional alerts that exist, press the UP ARROW or DOWN ARROW. All the "E" type errors are listed in **Section 7 - Troubleshooting**. Not all error sub-codes can be addressed by the user.

System Alert Icons:

Power was removed at indicated time

Power was restored at indicated time

More information on system alerts is provided in Section 4.12. A complete listing of system (E) and activity (W) alerts is provided in **Section 7 – Troubleshooting**.


Figure 5-25 – System Alerts Display

5.11 Needs Attention

The Needs Attention display enables you to view a list of panels in the instrument's database that have encountered a condition that requires operator attention. These conditions generally represent problems with the panels themselves, or with the information related to the panels. In many cases, the conditions can be corrected by operator action. See Figure 5-26.

When you view the Needs Attention panel list, the system provides the opportunity to resolve or ignore the condition that caused the panel to be placed in the list. If the panel was placed in the list due to missing or unresolved information (e.g., a Tie), the instrument provides the ability to add or edit the information to resolve the condition. If the panel was placed in the list due to a software, panel, or hardware error, the instrument provides the ability to delete the panel to resolve the error condition. Deleting a panel that is still testing causes that panel's protocol to be aborted. Only panels whose Needs Attention reason has not been ignored are shown in the display.

The Needs Attention display can be filtered to show one line per panel (filtered) or one line for each needs attention reason (unfiltered). A "toggle display" soft key enables you to switch the display between filtered (the default view) and unfiltered. When the display is filtered, only the highest priority reason is listed. Panels are shown in blue if additional needs attention conditions exist for that panel. When the display is unfiltered, each needs attention condition is listed separately. If you print a report, it prints either filtered or unfiltered, depending on the current view.

When you toggle from the unfiltered display to the filtered display, the highest priority Needs Attention reason code for the selected panel is highlighted. When you toggle back to the unfiltered display, the same reason code is highlighted if it still appears in the list (see below).

The Needs Attention display lists the first 30 panels with a Needs Attention flag, sorted by reason code (in the same order as the list below) and then by accession number. Panels without an accession number are listed first within each Needs Attention reason code. If more than 30 panels need attention, only the 30 panels with the highest priority are shown until some of the panels are resolved.

When you clear or ignore a Needs Attention reason code, the next highest priority reason code for that panel is highlighted. If there are no more reason codes for that panel, then the item at the top of the list is highlighted.


When you highlight a Needs Attention reason code and switch to the Panel Results display, when you return to the Needs Attention display, the same reason code is highlighted. If that reason is resolved, then the next highest priority reason code for that panel is highlighted. If there are no more reason codes for that panel, then the item at the top of the list is highlighted.

To access Needs Attention:

Press the "needs attention" soft key from the Main Status screen

Needs Attention soft keys:

"results"		"ignore"	
"delete panel"		"locate panel"	
"print"		"exit"	
"toggle display"		"clear NA"	

A panel may possess multiple Needs Attention reasons. If this is the case, the reasons are given in a prioritized list. The Needs Attention reasons are (in order of priority):

- **Test Aborted** – one of several incubator System Alerts occurred; or an ongoing panel that was in the instrument has not tested for more than one hour; or the instrument was turned off for more than one hour; or the instrument door has been opened for more than one hour; or system software did not execute testing algorithms for more than one hour; or a panel was moved to a different tier/instrument; or a media type was not specified for a Yeast ID panel
- **Can Not Identify Barcode** – internal barcode scanner could not read a panel sequence number in the same station where the instrument could determine a panel present condition (possible scratched or missing label); or invalid panel sequence number was scanned; or an unknown panel type was scanned
- **Can Not Read Panel Wells** – internal barcode scanner has read a sequence number in a station but the instrument cannot determine that a panel is present in that same station
- **Panel Lot Expired** – panel was logged in (or has a test start date) with a panel lot number that has already expired
- **Invalid AST Results** – at least one MIC field contains a value of "X"; excludes QC panels
- **Panel Missing** – internal barcode scanner reads a sequence number on a panel with an ongoing status during an inventory scan, but either the panel was removed before it completed its test protocol and hasn't obtained a result, or an internal scanner failure occurred and the sequence number can no longer be read
- **No Growth on Panel** – instrument did not detect growth in Growth Control well of panel; excludes QC panels and ID panels
- **Panel Lot Undefined** – Panel being logged in or placed in the instrument is from an undefined lot
- **Review QC Results** – the status of a completed QC panel is "REVIEW." Panel results should be reviewed to determine if an error occurred during panel preparation/handling (select status of REPEAT) or if QC test actually failed (select status of FAIL)
- **Missing Accession Number** – panel does not have any accession/isolate information; excludes QC panels

- **Missing Organism ID** – panel has no organism ID, or has unresolved tie IDs, and has no related panel with an organism ID (which includes related panels that have unresolved tie IDs or organism conflicts); excludes QC panels
- **Can Not Determine Organism ID** – panel has an Instrument Organism ID (or user-selected ID) of “No Identification” or it has no ID and its related panel has “No Identification” for its final ID, or ID panel concludes testing protocol without obtaining ID result; excludes QC panels
- **Invalid Organism ID** – panel has no defined organism ID and its related panels have a defined Organism ID of “Invalid Organism Received” (conferred to the panel via Auto Association); or panel receives download of an organism ID not in the Phoenix database; or panel is logged in or modified to “Invalid Organism Received” in the Organism ID field
- **Organism ID Conflict** – panel completes testing and has at least one related un-finalized panel that contains a different Final ID; excludes QC panels. Note that this code can be erroneously set as a result of conflicting system activities, such as when a panel completes testing and begins auto associating ID results as the system is modifying all the related panels with a user selected ID. If this occurs prior to the completion of the manually triggered modifications, the system can set this Needs Attention reason because IDs were in conflict due to two independent processes occurring. You can determine this by reviewing all panel ID results.
- **BDXpert Rule Flagged** – panel triggered one of the BDXpert rules defined, and the BDXpert rule is manually enabled in the Configuration menu, and the BDXpert rule triggered is pending; excludes QC panels, ID panels, and Combination panels with only the ID side enabled
- **Pending Too Long** – panel has not been scanned internally by an instrument (during an inventory scan) within 30 minutes of being entered in Panel Login

NEEDS ATTENTION				
Filtered on highest priority reason				
Reason	Sequence#	Accession#	Isolate	Station
Missing Organism ID	429530001405	1313	1	01/A05
Organism ID Conflict	429511111405	1313	1	01/A02

Figure 5-26 – Needs Attention Display

The following chart shows which soft keys are active for each reason code.

	"panel results"	"ignore"	"locate panel"	"delete panel"	"clear NA"
Test Aborted		*	*	*	
Can Not Identify Barcode		*	*		
Can Not Read Panel Wells			*	*	
Panel Lot Expired		*			*
Invalid AST Results		*			
Panel Missing		*		*	
No Growth on Panel		*	*	*	
Review QC Results	*	*	*	*	
Panel Lot Undefined					
Missing Accession Number	*		*		
Missing Organism ID	*				
Can Not Determine Organism ID	*	*			
Invalid Organism ID	*	*			
Organism ID Conflict	*				
BDXpert Rule Flagged	*				
Pending Too Long		*			

Table 5-3 – Needs Attention: Active Soft Keys

Needs Attention fields:

Note that some fields may be blank because panels are often placed in the Needs Attention list due to missing information.

Reason

Read-only field showing the Needs Attention reason (full list above)

Sequence #

Read-only field showing the panel barcode sequence number.

Accession #

Read-only field showing the accession number.

Isolate

Read-only field showing the isolate number.

Station

Read-only field showing the last known location of the panel, in the form II/Tst (where II is the instrument number, T is the Tier, and st is the station number).

Maintenance

6.1 General

The Phoenix system requires little maintenance from the user to provide reliable performance. Daily activities include checking the instrument temperature and printer paper supply. Weekly you should check the operation of the station status indicators (LEDs), the audible alarm, and the system alert indicator. Routine preventive maintenance consists of a monthly check of the recirculating air filters. All other procedures are on an "as needed" basis. Any maintenance or repair not described in this section should be performed by BD personnel only.

WARNING

ALL MAINTENANCE AND REPAIR OTHER THAN THE PROCEDURES DESCRIBED IN SECTION 6.2 – ROUTINE MAINTENANCE AND SECTION 6.3 – MODULE REPLACEMENT, MUST BE PERFORMED BY QUALIFIED SERVICE PERSONNEL. NON-COMPLIANCE WITH THIS WARNING MAY RESULT IN PERSONAL INJURY OR INSTRUMENT MALFUNCTION.

6.2 Routine Maintenance

6.2.1 Daily Maintenance

Each day the following maintenance procedures should be performed. The best time to perform maintenance is first thing in the morning, but it may be done any time you find convenient. The following items should be checked:

- 1 Check the printer's paper supply. If the paper supply is low or exhausted, replace the paper as explained in the manufacturer's operating instructions.
- 2 Record the temperature readout on the LCD display and on the temperature standard panel. Note that the temperature panel can be brought into view by selecting one of the instrument LED check functions on the Maintenance menu (see Section 6.2.2). This temperature should be $35^{\circ}\text{C} \pm 1.5^{\circ}\text{C}$.

6.2.2 Weekly Maintenance

Each week you should verify the operation of the instrument indicator lamps and audible alarm. (These items appear on the Daily Instrument Report and may be performed daily if desired.)

- 1 To perform lamp verification, access the Main Status screen (Figure 6-1):



Figure 6-1 – Main Status Screen

- 2 Press the "configuration" soft key:

- 3 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 4 From the main Configuration display, press the "maintenance" soft key.

- 5 Perform tests on the internal and external LEDs and audible alarm (see list below). To select an item, use the UP ARROW or DOWN ARROW key to highlight the item. Then press the "perform action" soft key (below). For the internal LED tests, after a moment, the "door unlocked" icon appears, and the corresponding tone sounds. Open the door (only required for the first item being checked ... the door may remain open after that until maintenance is complete) and verify the item being checked.

Each of the following items that appear in the maintenance menu should be checked:

Test internal green LEDs (causes green LEDs at the 20 accessible stations to illuminate)

Test internal red LEDs (causes red LEDs at the 20 accessible stations to illuminate)

Test internal amber LEDs (causes amber LEDs at the 20 accessible stations to illuminate)

Test System Alert Indicator – (causes system alert indicator to blink for five repetitions at one second on followed by one second off)

Test Alarm – (causes sample audible alarm tone to sound at the default alarm volume of 5)

If any item that is being verified (LEDs, audible alarm) fails to operate, contact your BD representative to schedule repair.

6.2.3 Periodic Maintenance

6.2.3.1 Air Filter Replacement

There are three air filters on the Phoenix instrument. The Electronics Bay Filter is located on the right side of the instrument above the I/O ports. The Source Bay Filter is located on the left side of the instrument just below the center. The Incubator Air filter is located on the right side of the center bay of the instrument, near the top.

Check and if necessary change these air filters every six months. If the instrument's environment is especially dusty, check the filters more frequently. These filters must remain clean and unobstructed; restricted air flow may cause the instrument interior to reach excessive temperatures, which can affect results and possibly cause hardware malfunctions or failures. It is particularly important to keep the incubator air filter clean.

Required Material

New air filters (see Appendix B for catalog numbers)

Medium phillips screwdriver (No. 2)

To Replace the Electronics Bay Air Filter (Figure 6–2)

- 1 The Electronics Bay Filter is located on the right side of the instrument above the I/O ports. The air filter is located behind a cover held in place by two phillips screws.
- 2 To remove the filter, use a No. 2 phillips screwdriver to loosen the two screws holding the filter cover in place. The screws are captive and require only one-quarter turn. They cannot be removed completely.
- 3 Flip the filter out and examine it for excessive dust and debris. Either replace the filter with a new one or return the existing filter back in the instrument.
- 4 Place the filter cover back in place. Tighten the two captive screws.

(Right side of instrument)


Figure 6-2 – Electronics Bay Filter Replacement

To Replace the Source Bay Air Filter (Figure 6-3)

- 1 The Source Bay Filter is located on the left side of the instrument just below the center. The air filter is located behind a cover held in place by two phillips screws.
- 2 To remove the filter, use a No. 2 phillips screwdriver to loosen the two screws holding the filter cover in place. The screws are captive and require only one-quarter turn. They cannot be removed completely.
- 3 Flip the filter out and examine it for excessive dust and debris. Either replace the filter with a new one or return the existing filter back in the instrument.
- 4 Place the filter cover back in place. Tighten the two captive screws.


Figure 6-3 – Source Bay Filter Replacement

To Replace the Incubator Air Filter (Figure 6-4)

- 1 The Incubator Air filter is located on the right side of the center bay of the instrument, near the top. The air filter is located behind a cover held in place by two phillips screws.
- 2 To remove the filter, use a No. 2 phillips screwdriver to loosen the two screws holding the filter cover in place. The screws are captive and require only one-quarter turn. They cannot be removed completely.
- 3 Flip the filter out and examine it for excessive dust and debris. Either replace the filter with a new one or return the existing filter back in the instrument.

CAUTION

If you return the existing filter, insure that you replace the filter in the same orientation that it was removed (i.e., with arrow stamped on the side of the filter pointing in the same direction).

- 4 Place the filter cover back in place. Tighten the two captive screws.


Figure 6-4 – Incubator Filter Replacement

6.2.3.2 Cleaning the Barcode Scanner Window

There are no user-serviceable parts in the barcode scanners. After a period of use, you may find that a barcode scanner does not seem to scan barcode labels as easily as before. If this happens, you can try to solve the problem by cleaning the scanner's window. To clean the window, use a lint-free, nonabrasive cloth dampened with distilled water. Dry the window with a dry lint-free nonabrasive cloth.

6.2.3.3 Cleaning/Decontamination

A situation requiring biological decontamination of one or more panel locations can occur if a panel should leak while in the instrument. The priority in this situation is to first limit the extent of the contamination and then to decontaminate the panel location(s) and other accessible instrument areas receiving the spill. If the spill extends into regions of the carousel not accessible for topical decontamination, contact BD Technical Service (USA) for further instructions (800-638-8663).

To Decontaminate Carousel Panel Locations

The solution recommended to clean the affected surfaces should be at least a 10 percent household bleach solution. All surfaces must be thoroughly washed with the freshly prepared bleach solution, so that the surfaces are "glistening wet." If you are not sure of the extent of the contamination, thoroughly wash the exposed portions of the carousel and cabinet with the freshly prepared bleach solution.

WARNING

ALL PORTIONS OF THE BODY THAT COULD POSSIBLY COME INTO CONTACT WITH THE AFFECTED INSTRUMENT SURFACES MUST BE COMPLETELY COVERED BEFORE BEGINNING THE DECONTAMINATION PROCESS.

- 1 Wear gloves and a gown, completely covering any body surfaces that could possibly come into contact with the affected instrument surfaces.
- 2 Turn off power to the instrument. Unplug the instrument power cord before proceeding.
- 3 Completely absorb the contaminated spill (gauze pads are most effective).
- 4 Apply the bleach solution to the affected surfaces, so that the surfaces are "glistening wet." Let stand for approximately 15 minutes.
- 5 Absorb the applied solution with gauze pads or paper towels.
- 6 Dampen a clean cloth with water. Wipe down the decontaminated surfaces.
- 7 Thoroughly dry all wet surfaces.
- 8 Discard all cleanup materials as biohazardous waste.

6.2.4 Maintenance Menu

The Maintenance menu, accessible from the Configuration menu, provides several functions for performing instrument maintenance. A number of these functions, noted below, are for BD use only. Others are for user weekly and "as needed" maintenance. The following functions appear under the Maintenance menu:

- Test Internal Green LEDs (Section 6.2.2)
- Test Internal Red LEDs (Section 6.2.2)
- Test Internal Amber LEDs (Section 6.2.2)
- Test System Alert Indicator (Section 6.2.2)
- Test Alarm (Section 6.2.2)
- Install/Upgrade (Section 2.5)
- Save System Data to Disk (Section 7.3)
- Save Syslog to Disk (Section 7.4)
- Access Normalizer Panel (for BD use only)
- Save User Data to Disk (Section 6.2.4.1)
- Restore User Data from Disk (Section 6.2.4.2)
- Save LIS Codes to Disk (Section 6.2.4.3)
- Restore LIS Codes from Disk (Section 6.2.4.4)
- Save Panel Lot Definitions (Section 6.2.4.5)
- Restore Panel Lot Definitions (Section 6.2.4.6)
- Save Panel Configuration (Section 6.2.4.7)
- Install Panel Configuration (Section 6.2.4.8)
- Phoenix Update Disk (Section 6.2.4.9)
- Modify Panel Usage (Section 6.2.4.10)
- Write Lightsrc to Diskette (Section 6.2.4.11)
- Read Lightsrc from Diskette (Section 6.2.4.12)
- View Syslog (Section 6.2.4.13)
- Save Syslog to Network (Section 6.2.4.14)
- Display UV and RGB Status (Section 6.2.4.15)
- Force Front UV Bulb Adjustment (Section 6.2.4.16)
- Force Rear UV Bulb Adjustment (Section 6.2.4.17)
- Force RGB LED Adjustment (Section 6.2.4.18)
- Front Panel Reader Setup (Section 6.2.4.19)

6.2.4.1 Save User Data to Disk

Save User Data to Disk enables you to backup the configuration parameters you have set to floppy disk. This data includes: Custom Interpretation Rule Set Configuration, and BDXpert Rules Configuration. The information saved is for the current instrument only. However, information saved at one instrument can be restored at another instrument. This enables you to have consistent configurations among all the instruments. To save data to disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE USER DATA TO DISK.*"
- 5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.
- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to save the data to diskette. Press the "confirm" soft key to proceed with the save or "cancel" to cancel the operation.

OK

or

- 8 An hourglass icon appears while the information is being written to the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

- 10 Write protect the disk.

6.2.4.2 Restore User Data from Disk

Restore User Data from Disk enables you to restore the configuration parameters you have copied to floppy disk back to the instrument. This data includes: Custom Interpretation Rule Set Configuration, and BDxpert Rules Configuration. To restore, the instrument door must be closed. To restore data from disk, follow the steps below:

NOTE

This operation should not be performed while there are ongoing panels.

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "RESTORE USER DATA FROM DISK."
- 5 Insert the saved user data floppy disk in the disk drive.
- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to restore the data from diskette, and advises you that the instrument will reboot when the restore is complete. Press the "confirm" soft key to proceed with the restore or "cancel" to cancel the operation.

OK

or

- 8 The system verifies the dates of the backed up files on the floppy disk. If any files are older than the ones currently on the instrument, you will be asked to confirm that you want to replace the more current files on the instrument with the older ones on the floppy disk. You can choose to replace none, one, several, or all of the newer files.
- 9 An hourglass icon appears while the information is being copied from the disk, as well as the message, "Restoring Data."

- 10 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. The instrument then automatically reboots.

- 11 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button. **Note that the new files may take at least 15 minutes to load into instrument memory.**

6.2.4.3 Save LIS Codes to Disk

Save LIS Codes to Disk enables you to save all the Organism and Antimicrobial LIS codes you have defined to floppy disk. This produces a text file that you can edit on a PC, which might be quicker for some users who have many edits to perform. Then, the edited codes can be restored (Section 6.2.4.4) back to the Phoenix instrument. The function also enables you to copy codes from one instrument to another. Save LIS Codes to Disk appears on the Maintenance menu only if LIS Communications is enabled.

Entries in the text field consist of:

Identifier (ORG, QC_ORG, DRUG) | BD Code | Short Name/Drug | LIS Code

Codes cannot be added or deleted, and **only** the LIS Code portion can be modified. If another field is changed, the instrument will not restore the codes.

To save codes to disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.


- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE LIS CODES TO DISK."
- 5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.
- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to save the data to diskette. Press the "confirm" soft key to proceed with the save or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the information is being written to the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.4 Restore LIS Codes from Disk

Restore LIS Codes from Disk enables you to restore the Organism and/or Antimicrobial codes you previously copied to floppy disk back to the instrument. The restore operation completely overwrites the existing Organism/Antimicrobial LIS Code database. Codes will not be restored if any field other than LIS Code was modified, or if LIS Code was entered in an incorrect format. If this happens, an error log is written on the floppy disk, and you can review this file to see what caused the error. Restore LIS Codes from Disk appears on the Maintenance menu only if LIS Communications is enabled.

To restore codes from disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "RESTORE LIS CODES FROM DISK."

- 5 Insert the saved LIS codes floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to restore the data from diskette. Press the "confirm" soft key to proceed with the restore or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the information is being copied from the disk, as well as the message "Restoring Data."

- 9 When the activity is complete, the message, "Complete" flashes on the screen, the "activity completed" tone sounds, and the hourglass icon disappears.
- 10 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.5 Save Panel Lot Definitions

Save Panel Lot Definitions enables you to save the panel lots present in the instrument to a floppy disk. This enables you to transfer these lot definitions (and QC panel results) to other instruments so that the records can be viewed/used there. **Save Panel Lot Definitions** saves the defined panel lots and related data (Sequence #s, Expirations Dates, etc.), and QC results for any strains tested in those lots. Save Panel Lot Definitions appears on the Maintenance menu only if QC Lot Support is enabled.

To save panel lot definitions to disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE PANEL LOT DEFINITIONS."
- 5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.
- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to save the data to diskette. Press the "confirm" soft key to proceed with the save or "cancel" to cancel the operation.

OK

or

- 8 An hourglass icon appears while the information is being written to the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.6 Restore Panel Lot Definitions

Restore Panel Lot Definitions enables you to restore Panel Lot Definitions saved at one instrument to another. This makes Panel Lot definitions *and* QC panel results transferrable in labs that use multiple instruments, so that you only have to define a lot once using the box (carton) label. Restore Panel Lot Definitions appears on the Maintenance menu only if QC Lot Support is enabled.

To restore panel lot definitions from disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "RESTORE PANEL LOT DEFINITIONS."
- 5 Insert the saved definitions floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to restore the data from diskette. Press the "confirm" soft key to proceed with the restore or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the information is being copied from the disk, as well as the message "Restoring Data."

- 9 If you receive a message about a conflict between the restore disk and the instrument's panel lot information, contact BD for assistance.

Otherwise, when the activity is complete, the message, "Complete" flashes on the screen, the "activity completed" tone sounds, and the hourglass icon disappears.

- 10 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.7 Save Panel Configuration

Save Panel Configuration enables you to backup the instrument's panel configurations to floppy disk. The information saved is for the current instrument only. However, information saved at one instrument can be installed at another instrument. This enables you to have consistent panel configurations among all the instruments. To save panel configurations to disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE PANEL CONFIGURATION."
- 5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to save the data to diskette. Press the "confirm" soft key to proceed with the save or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the information is being written to the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk drive indicator is off, then remove the floppy disk by pressing the disk eject button.
- 10 Write protect the disk.

6.2.4.8 Install Panel Configuration

Install Panel Configuration enables you to update the instrument's panel configurations. To perform this operation, the instrument door must be closed.

To install panel configurations, follow the steps below:

NOTE

This operation should not be performed while there are ongoing panels.

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "INSTALL PANEL CONFIGURATION."
- 5 Insert the (first, if there is more than one) panel configuration floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to install the data from diskette, and advises you that the instrument will reboot when the installation is complete. Press the "confirm" soft key to proceed with the installation or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the information is being copied from the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. The instrument then automatically reboots.
- 10 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.9 Phoenix Update Disk

Phoenix Update Disk (PUD) enables you to update numerous data files in the instrument, such as antimicrobial breakpoints and rules, QC data, drugs, organisms, etc. To install the Phoenix Update Disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "PHOENIX UPDATE DISK."
- 5 Insert the (first, if there is more than one) update disk in the disk drive.
- 6 Press the "perform action" soft key.


- 7 The system asks you to confirm that you wish to install the data from diskette, and advises you that the instrument will reboot when the installation is complete. Press the "confirm" soft key to proceed with the installation or "cancel" to cancel the operation.

OK

or

- 8 The system verifies the version and configuration of the files on the floppy disk. If any files are older than the ones currently on the instrument, the instrument issues an alert (E34) and does not perform the update.
- 9 An hourglass icon appears while the information is being copied from the disk.

- 10 If additional disks are required, the system displays the following message: "Please insert disk # n of Phoenix update set, then press OK to continue."

Remove the disk and insert Phoenix update floppy disk # n in the disk drive.

- 11 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. The instrument then automatically reboots.
- 12 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

6.2.4.10 Modify Panel Usage

The Modify Panel Usage operation enables you to disable either the ID or the AST side of a combination panel that is not Removable (provided both sides are currently enabled). To modify panel usage:


- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the “perform action” soft key.

- 3 From the main Configuration display, press the “maintenance” soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, “MODIFY PANEL USAGE.”
- 5 The following display appears:

MODIFY PANEL USAGE

Sequence #

Accession #

Isolate #

☐

ID

☒

AST

☒

- 6 Enter or scan the panel barcode sequence number of the panel whose usage you wish to modify. The system automatically completes the **Accession #** and **Isolate #** fields, which are read-only.
- 7 Use the "tab" key to advance to either the **ID** or **AST** Checkbox (depending on which side of the panel you want to disable). Press the "select" soft key or Spacebar to "un-check" the checkbox.
- 8 Press the "save" soft key to save the panel modification.

Note the following conditions for panel usage modification:

- ◆ No information for related panels is modified when panel usage is changed.
- ◆ If there is no user-entered Final ID and the ID side is disabled, the Instrument ID, Biochemical Results, Confidence Values, SIR Values, and ID Special Messages are removed from the record.
- ◆ If there IS a user-entered Final ID and the ID side is disabled, this ID is retained, as are SIR Values. However, Instrument ID, Biochemical Results, Confidence Values, and ID Special Messages are removed from the record.

6.2.4.11 Write Light Source to Disk

The Write Light Source to Disk operation enables you to create a floppy disk containing the current values needed by the instrument to establish setpoints for source LEDs. **The operation should only be performed when advised by your local BD representative.**

- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.


- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Insert a blank, formatted, write-enabled floppy disk in the disk drive. Be sure the disk is labeled "Light Source Disk" and includes the instrument serial number.
- 5 Use the UP ARROW or DOWN ARROW to highlight the selection, "WRITE LIGHTSRC TO DISKETTE." The instrument verifies that a proper disk is in the drive.
- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to save the data to diskette. Press the "confirm" soft key to proceed with the save or "cancel" to cancel the operation.

- 8 An hourglass icon appears while the Light Source information is being written to the disk.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk drive indicator is off, then remove the floppy disk by pressing the disk eject button.
- 10 Write protect the disk.

6.2.4.12 Read Light Source from Diskette

Read Light Source from Disk enables you to install new light source files that have been saved on floppy disk (either the one provided with the system at installation or one that BD directed you to create). The operation can be performed only when the instrument is idle, there are no ongoing panels, and the instrument has completed an inventory scan. The same normalizer panels must be present in the instrument as when the light source disk was made.

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "READ LIGHTSRC FROM DISKETTE."
- 5 Insert the light source floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 The system asks you to confirm that you wish to write the data from diskette to NVRAM. Press the "confirm" soft key to proceed with the installation or "cancel" to cancel the operation.

- 8 The system verifies the dates of the files on the floppy disk. If any files are older than the ones currently on the instrument, you will be asked to confirm that you want to replace the more current files on the instrument with the older ones on the floppy disk. You can choose to replace none, one, several, or all of the newer files.
- 9 An hourglass icon appears while the information is being copied from the disk.

- 10 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears.
- 11 Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.


6.2.4.13 View Syslog

View Syslog enables you to view a list of messages that the LIS communications function has generated. Note that the display will contain no messages if the BD EpiCenter system is connected. The messages represent status messages that have occurred during communications, such as query messages, log entries, interface messages, etc. This list can be filtered by date, printed (a potentially large report that prints in landscape orientation), and can be searched to find particular words or phrases (with repeat find of the term).

Refer to Appendix J – Event Log Messages for a complete listing of messages and corrective actions.

View Syslog soft keys:

"print"

"find"

"repeat data"

"select filter"

"exit"

To view the system log:

- 1 Press the "configuration" soft key.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "VIEW SYSLOG" and press the "perform action" soft key. The system log then appears.
- 5 All available LIS messages appear in the window, with the oldest messages at the top. If more messages exist than can fit in the window, use the DOWN ARROW to view messages below those displayed. If the message is wider than the window size, press  on the keyboard to see the rest of the message. Press  to return the display to the original position.

To search the system log:

- 1 Press the "find" soft key.

- 2 The Find display appears.
- 3 In the Find field, you can enter up to 37 characters (not case sensitive) for which you want the system to search.
- 4 The Direction checkboxes enable you to specify whether the search should proceed upward or downward (the default). Press the "select" soft key to check or uncheck the desired direction. Only one direction can be selected.
- 5 Press the "perform action" soft key to begin the search. The display returns to the system log, with the first message containing the term at the top of the window and highlighted. If the term is not found, no messages are highlighted.
- 6 To find the next occurrence of the term, press the "repeat data" key. If the term appears again, it is highlighted in the system log. If the term does not appear again, the same message remains highlighted and the "repeat data" key disappears.

To filter the system log:

- 1 Press the “select filter” soft key.

- 2 You may enter a Start Date/Time, an End Date/Time, both, or neither.

If you enter only a start date, the filter returns all messages from that date forward to the current time.

If you enter only an end date, the filter returns all messages from that date back to the oldest message in the log.

If you enter both dates, the filter returns all messages between the two dates.

If you enter neither date, all messages in the log are displayed (same as when the log is first accessed).

- 3 The Filter Selection display appears.

MAINTENANCE							
Dates ☒ Start Date 07 / 30 / 01 09 : 20 ☐ End Date							
↑		☞	○ ● ○				☛

- 4 To enter a Start Date/Time or End Date/Time, highlight the desired checkbox and press the "select" soft key. Date and time fields appear when the box is checked. Move the cursor to the day, month, year, hour, or minute field, and press the UP ARROW or DOWN ARROW to increase or decrease the displayed value.
- 5 When the desired dates and times are entered, press the "perform action" soft key to execute the filter.

6.2.4.14 Save Syslog to Network

Save Syslog to Network should only be used when advised by your local BD representative.

The **Save Syslog to Network** option appears only when the instrument is connected to a BD EpiCenter system. This option enables you to save the event log to an EpiCenter system.

Under certain circumstances, a BD representative may advise you to save (write) the syslog to the EpiCenter system. These circumstances include some error conditions and system malfunctions. The Save Syslog function copies the system event log, which contains logged system messages about various system, instrument, and communications events. The instrument door must be closed, and the instrument must be idle to save the syslog.

- 1 Press the "configuration" soft key.

- 2 The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE SYSLOG TO NETWORK."
- 5 Insure that the data collection application is launched on the EpiCenter computer.
- 6 Press the "perform action" soft key.

- 7 An hourglass icon appears while the information is being sent.

- 8 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears, and the event log file is saved to the directory in which the data collection application resides on the EpiCenter computer.

6.2.4.15 Display UV and RGB Status

The Display UV and RGB Status operation enables you to view information about the UV and RGB source bulbs. **The operation should only be performed when advised by your local BD representative.**


- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Use the UP ARROW or DOWN ARROW to highlight the selection, "DISPLAY UV AND RGB STATUS."

- 5 Press the "perform action" soft key.

- 6 The instrument displays the following information on the left side of the screen:

For Front and Rear UV Bulbs:

In Service (YES/NO)
 Telsum (value)
 Cycles (value)
 Time (Min) (value)
 Available (YES/NO)
 Intensity (%)
 Adjustment (-----/PEND)

For RGB Source LEDs:

Adjustment (-----/PEND)

- 7 Press the "exit" soft key to return to the Maintenance menu.

6.2.4.16 Force Front UV Bulb Adjustment

The Force Front UV Bulb Adjustment orders the instrument to adjust the front UV bulb. **The operation should only be performed when advised by your local BD representative.** If the front UV bulb is currently in use, a message appears informing you that the adjustment will not occur until there are no ongoing panels in the instrument. If the rear UV bulb is currently in use, the front UV bulb will become the "in service" bulb and the adjustment will not occur until there are no ongoing panels in the instrument.

- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.


- 2 The password window appears. Enter the appropriate password and press the "perform


action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Use the UP ARROW or DOWN ARROW to highlight the selection, "FORCE FRONT UV BULB ADJUSTMENT."
- 5 Press the "perform action" soft key.

- 6 If the front UV bulb is currently in service, a message appears informing you that the adjustment remains pending until there are no ongoing panels in the instrument.

If the rear UV bulb is currently in service, a message appears informing you that the front UV bulb will become the "in service" bulb. The adjustment remains pending until there are no ongoing panels in the instrument.

To continue the UV adjustment process, press the "confirm" soft key. To cancel the process, press the "cancel" soft key.

OK

or

- 7 The system waits until there are no ongoing panels and then performs the UV bulb adjustment.
- 8 An hourglass icon appears while the adjustment is being performed.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears.
- 10 Press the "exit" soft key to return to the Maintenance menu.

6.2.4.17 Force Rear UV Bulb Adjustment

The Force Rear UV Bulb Adjustment orders the instrument to adjust the rear UV bulb. **The operation should only be performed when advised by your local BD representative.** If the rear UV bulb is currently in use, a message appears informing you that the adjustment will not occur until there are no ongoing panels in the instrument. If the front UV bulb is currently in use, the rear UV bulb will become the "in service" bulb and the adjustment will not occur until there are no ongoing panels in the instrument.

- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.

- 2 The password window appears. Enter the appropriate password and press the "perform

action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Use the UP ARROW or DOWN ARROW to highlight the selection, "FORCE REAR UV BULB ADJUSTMENT."
- 5 Press the "perform action" soft key.

- 6 If the rear UV bulb is currently in service, a message appears informing you that the adjustment remains pending until there are no ongoing panels in the instrument.

If the front UV bulb is currently in service, a message appears informing you that the rear UV bulb will become the "in service" bulb. The adjustment remains pending until there are no ongoing panels in the instrument.

To continue the UV adjustment process, press the "confirm" soft key. To cancel the process, press the "cancel" soft key.

OK

or

- 7 The system waits until there are no ongoing panels and then performs the UV bulb adjustment.
- 8 An hourglass icon appears while the adjustment is being performed.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears.
- 10 Press the "exit" soft key to return to the Maintenance menu.

6.2.4.18Force RGB LED Adjustment

The Force RGB LED Adjustment orders the instrument to adjust the Red, Green, and Blue source LEDs. **The operation should only be performed when advised by your local BD representative.** When you select the operation, a message appears informing you that the adjustment will not occur until there are no ongoing panels in the instrument.


- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Use the UP ARROW or DOWN ARROW to highlight the selection, "FORCE RGB LED ADJUSTMENT."

- 5 Press the "perform action" soft key.

- 6 A message appears informing you that the adjustment remains pending until there are no ongoing panels in the instrument. To continue the adjustment process, press the "confirm" soft key. To cancel the process, press the "cancel" soft key.

OK

or

- 7 The system waits until there are no ongoing panels and then performs the RGB LED adjustment.
- 8 An hourglass icon appears while the adjustment is being performed.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears.
- 10 Press the "exit" soft key to return to the Maintenance menu.

6.2.4.19 Front Panel Reader Setup

The Front Panel Reader Setup enables you to view and/or change how the instrument's built-in barcode scanner is treated when the instrument reboots. **The operation should only be performed when advised by your local BD representative.** The 2 choices offered are to 1) reset the scanner to factory default settings, or 2) retain the current settings when the instrument reboots.

- 1 From the Main Status Screen, press the "configuration" soft key. The main configuration screen appears.

- 2 The password window appears. Enter the appropriate password and press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Use the UP ARROW or DOWN ARROW to highlight the selection, "FRONT PANEL READER SETUP."

- 5 Press the “perform action” soft key.

- 6 The current Startup Status is shown at the left of the display.
- 7 Press the “perform action” soft key to change the displayed setting. Set the Startup Status value to “Initialize” to reset the scanner to factory default settings when the Phoenix instrument reboots. Set the Startup Status value to “Retain” to allow the scanner to retain its current settings when the instrument reboots.
- 8 Press the “exit” soft key to return to the Configuration menu.

6.3 Module Replacement

6.3.1 General

The Phoenix instrument has been designed and tested for troublefree performance. In the event of malfunction, you should contact BD for service under the terms of your existing service contract or warranty. Only the external barcode scanner is user-replaceable.

Replacement modules may be swapped for failed modules which are then returned to BD. Credit is then applied towards the replacement module. Only replacement parts supplied by BD should be used in the procedures described in this section.

6.3.2 Barcode Scanner Replacement

The barcode scanner is connected to the instrument at the I/O panel on the right side of the instrument. Refer to Figure 6-5.

Required Materials

1/8 in. flat screwdriver

To Replace the Barcode Scanner

- 1 On the I/O panel on the right side of the instrument, locate the cable connection for the barcode scanner. This cable is held in place by two small screws. Unscrew the screws and unplug the cable.
- 2 Note that the cable connector is roughly D-shaped, with the top of the connector being wider than the bottom part. Plug in the new scanner cable, being sure to align the cable connector to the plug in the correct orientation. Tighten the screws.
- 3 To verify proper operation of the new scanner, go to the Panel Login screen, advance the cursor to the **Sequence #** field, and scan a panel number barcode. Verify that the correct number is entered in the field. Press the "exit" soft key to exit the display without saving your entry.


Figure 6-5 – Barcode Scanner Replacement

6.3.3 Thermometer Removal

If you find that the fluid in the thermometer of the temperature standard panel has separated, follow the procedure below to remove the thermometer and replace it.

To Remove/Replace the Thermometer

- 1 Remove the temperature standard panel from the Phoenix instrument.
- 2 Using the small access hole at the bottom of the panel, gently push the thermometer upward through the large slotted opening at the top of the temperature panel.
- 3 Using your fingers pull the thermometer completely out of the panel through the slotted opening at the top of the temperature panel.
- 4 Reunite the separated fluid column per the instructions below.
- 5 Install the thermometer in the reverse order described above.

6.3.4 Reuniting Separated Liquid in the Thermometer

If you find that the fluid in the thermometer of the temperature standard panel has separated, follow the procedure below to reunite the liquid.

To Reunite Separated Liquid

- 1 Should there be a separation in the capillary or in the expansion chamber at the top of the thermometer, heat the bulb of the thermometer in a hot liquid which exceeds the range of the thermometer until the separation and main liquid column enter the expansion chamber and unite with each other.
- 2 Quickly remove the thermometer from the liquid, so that the liquid does not completely fill the expansion chamber which could possibly harm the thermometer.
- 3 Check the thermometer against a certified, traceable thermometer or an ice bath to assure the thermometer is reading correctly.

Troubleshooting

7.1 General

7.1.1 Instrument Service

If your BD Phoenix 100 instrument malfunctions or operates unusually in any way, you may initially attempt to solve the problem by following the recommendations in this section. All other servicing attempts will terminate the responsibility of the manufacturer under the terms of the warranty.

If you cannot repair a system malfunction, contact your local BD representative (contact numbers are listed in Appendix D).

This section discusses error messages, which appear when the system has encountered a known problem. These messages are listed in numerical order, along with possible causes of the message and corrective actions.

7.2 Error/Alert Messages

CAUTION

When the system notifies you of alerts and errors, you should immediately respond to the condition.

When the system encounters an alert or error condition, the error code (EXX or WXXX, where XX or XXX is a number) is either displayed on the screen or written into the system alert list. The error code is an abbreviation for the conditions described in the listing below.

Different types of alerts and errors behave in different ways. There are three basic types of alert condition:

- ① Self-Clearing – these alerts are removed from the System Alert Screen after they have been displayed and the screen has been exited. No other intervention is required.
- ② Persistent – these alerts remain in the System Alert Screen (even after displaying the alert code) until the instrument determines that the error causing the alert has been corrected. Correction can be accomplished with or without user intervention, depending on the alert.
- ③ Auto Clear – these alerts are removed from the System Alert Screen as soon as the condition causing the error has been cleared. There is no requirement to display the alert code on the System Alerts List Screen. If the alert is still present when reviewing the Alert list, it will be displayed.

"W" error codes are displayed on the screen when they occur. (They also cause the Activity Error tone to sound [sequence of short high beep and short low beep repeated four times].) These are activity (or "workflow") types of errors. In most cases, this means that some action you have performed was not what the system expected, but you can usually perform the correct action, as recommended below, without exiting the current operation. These activity errors are flagged by the Activity Error icon:

"W" error codes are grouped as follows:

W1XX	Problem with the Sequence Number
W2XX	Problem with the Accession, Isolate Number, or Media type
W3XX	Problem with a QC panel only
W4XX	Action to get into the instrument is not allowed
W5XX	Configuration/Maintenance screen activity not allowed (as long as it is not an instrument action)
W6XX	Screen activity is not allowed
W7XX	LIS errors
W8XX	QC Lot Support errors

System alerts, which comprise all "E" type error codes, are reported in the system alert list. These errors cause the Alert tone (medium beep on for one second, off for 3 seconds, repeating) to sound (if it is enabled). Also the System Alert icon appears on the Main Status screen and the System Alert indicator flashes. The errors must be reviewed to clear the system alert condition. The system alert list can be viewed from the Main Status screen by pressing the "system alert" soft key:

The error codes are listed in numerical order. Error sub-codes are 8-digit numeric codes that appear below the EXX readout in the system alert list. The sub-codes indicate specific conditions detected. Many sub-codes are listed in the "Possible Causes" and "Corrective Actions" sections below.

CAUTION

If any error sub-codes other than those listed here appear, note the sub-code and contact BD for assistance.

If the recommended corrective actions do not solve the problem, contact BD.

E01 Incubator temperature alert

Temperature readings are taken every 5 seconds for 10 minutes and these readings are averaged.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

00000001– Average temperature too high (> 36.5° C)

00000002 – Average temperature too low (< 33° C)

00000010 – Average temperature too high for more than one hour (> 36.5° C)

00000020 – Average temperature too low for more than one hour (< 33° C)

00000040 – Absolute high temperature (> 38.5° C)

Note that if the instrument temperature reaches 39° C ± 0.5° C, the instrument disables the heater.

For all sub-codes, causes include: room temperature is not within recommended range, or other environmental specification is not being met (such as instrument sitting in direct sunlight or too close to HVAC air register). Air filters are dirty, restricting fresh air intake.

CORRECTIVE ACTION(S)

- ☐ Make sure environmental specifications are met (see Section 2 – Installation and Setup).
- ☐ Replace air filters.

E02 Lightbay temperature alert

The instrument has detected either that the temperature is too high or too low on any of the tiers, or communications with the lightbay sensor is terminated. The instrument continues testing on schedule.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Temperature too High (Light Source: > 50° C; Tier Micro: > 60° C)
Tier A (Applicable subcodes: 00000001, 00000008)
Tier B (Applicable subcodes: 00000040, 00000200)
Tier C (Applicable subcodes: 00001000, 00008000)
Tier D (Applicable subcodes: 00040000, 00200000)
- Temperature too Low (< 18° C):
Tier A (Applicable subcodes: 00000002, 00000010)
Tier B (Applicable subcodes: 00000080, 00000400)
Tier C (Applicable subcodes: 00002000, 00010000)
Tier D (Applicable subcodes: 00080000, 00400000)
- Causes include: room temperature is not within recommended range, or other environmental specification is not being met (such as instrument sitting in direct sunlight or too close to HVAC air register). Air filter is dirty, restricting fresh air intake.

CORRECTIVE ACTION(S)

- ☐ Make sure environmental specifications are met (see Section 2 – Installation and Setup).
- ☐ Replace air filters.

I 13 Ambient temperature alert

Communication with the ambient sensor is terminated. The instrument continues testing on schedule.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Failure of ambient temperature sensor

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E05 Carousel alert

The instrument reports this alert if one of the following occurs: the reported RPM is below spec, certain flag readings are incorrect, or the carousel is jammed or stalled. It is detected during any carousel rotation. Any inventory scan or panel test in progress is aborted and the instrument ignores any data received from the test or scan.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- 00000001 – RPM is below spec
- 00000040 – Carousel is stalled
- 00000080 – Carousel is jammed

Listed codes: something is impeding the motion of the carousel.

CORRECTIVE ACTION(S)

- ☐ Open the door, look for and remove any obstructions such as a panel that is ajar or a panel closure that is not seated. **Do not manually rotate the carousel.** Close the door. If message recurs, contact your local BD representative.

E06 Tier alert

This is the general alert condition for problems related to specific tiers. Such problems include: flag readings are incorrect, optical errors, timeout conditions, and normalizer problems. The alert is reported in the following format (typical sub-codes are shown):

E06	00000001
	00000000
	00000100
	00000280

The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

Tier alerts are detected during a test cycle or inventory scan. The instrument ignores any data received from the test or scan for the tier that has sustained the error. Data received from good tiers is retained. All stations in the bad tier are automatically blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- 00000020 – Rotor timeout
- 00002000 – Tier out of service – the instrument cannot communicate with the tier

CORRECTIVE ACTION(S)

- ☐ 00000020 – Check for/remove object impeding motion of carousel
- ☐ 00002000 – Reboot instrument (turn power off, wait 10 seconds, turn power back on)

E07 Power supply alert

One of the power supplies was out of specification. The instrument ignores any data received from the test.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Temporary electrical anomaly.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E08 Tier power supply alert

One of the tier power supplies was out of specification. The instrument ignores any data received from the test for panels in that tier. The alert is reported in the following format (typical sub-codes are shown):

```
E17  00000001
      00000002
      00000004
      0000000A
```

The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

Tier alerts are detected during a test cycle. The instrument ignores any data received from the test for the tier that has sustained the error. Data received from good tiers is retained. All stations in the bad tier are automatically blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Temporary electrical anomaly.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E09 Test aborted

Panel testing has not occurred for more than one hour.

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- 20000000 – Causes include: instrument was off or door was open for more than an hour, system clock was set ahead one hour or more, test cycle has not occurred in more than an hour, etc.

CORRECTIVE ACTION(S)

- ☐ All ongoing panels are set to Needs Attention, and their status is set to Removable. Some results may be incomplete, and all affected panel results should be reviewed.

E10 System database corruption

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- 00000002, 00000008, 00000020, 00000040, 00000100, 00000200, 00000400 – Database checksum test failed or data in one of the databases (see below) was corrupted and replaced by a default database.

CORRECTIVE ACTION(S)

- ☐ Write data to disk (see Section 7.3) and call BD.
- ☐ 00000002 – Check your settings in the Configuration screens and reset them to your preferences (see Section 2.3.1).
- ☐ 00000008 – Check your settings in BDxpert Rule Configuration and reset them to your preferences (see Section 2.3.4).
- ☐ 00000020 – Check your breakpoint settings in Custom Interpretation Rule Set Configuration and reset them to your preferences (see Section 2.3.3).
- ☐ 00000040 – Check your settings in Organism and Antimicrobial Configuration and reset them to your preferences (see Sections 2.3.2.1 and 2.3.2.2).
- ☐ 00000100 – Perform a Read Light Source from Diskette operation as described in Section 6.2.4.12.
- ☐ 00000200 – Perform a Read Light Source from Diskette operation as described in Section 6.2.4.12.
- ☐ 00000400 – Informational message indicating that LIS orders have been deleted from the database. This can occur when there is a checksum error in the LIS orders area of the database, or the internal data structures have been changed during a software update. You should resend any orders that were deleted as a result of this alert.

E11 Printer error

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

20000000

- Printer paper is jammed or exhausted.
- Printer cable is disconnected during printing.
- Printer power is turned off during printing.
- Printer was taken offline during printing.

CORRECTIVE ACTION(S)

- ☐ Check paper and clear jam or add paper if necessary.
- ☐ Check printer cables (power and communications), reattach if necessary.
- ☐ Turn printer power on.
- ☐ Place printer online.
- ☐ Report will resume printing automatically when error condition is corrected.

E12 Source monitor error

This is the alert condition for problems related to testing LEDs in the tiers. The alert is reported in the following format (typical sub-codes are shown):

E12	00000100
	00000000
	00000000
	00000000

The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

Source monitor errors are detected during a test cycle. The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Optical problem with source LEDs.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E13 Power failure

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- 20000000 – Power was removed from instrument.

CORRECTIVE ACTION(S)

- ☐ Message is informational. If multiple power failures have occurred, only the latest one is reported in the alert list. Note the power failure and restore times in your instrument log. Note that power fail events are not recognized until the instrument user interface has successfully loaded.

E14 CCD underrun

Scanning of a panel stopped prematurely. The alert is reported in the following format:

```
E14      00000000
          00000100
          00000000
          00000000
```

The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

CCD underrun errors are detected during a test cycle. The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Problem in the tier.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E15 UV bulb switched

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

Either:

- 00000001 The rear ultraviolet bulb
- 00000002 The front ultraviolet bulb

was switched.

CORRECTIVE ACTION(S)

- ☐ Message is informational. The switch happened as a result of an automatic light source adjustment.

E16 485 master communication error

There was a communications error relating to testing. The alert is reported in the following format:

E16	00000001
	00000000
	00000000
	00000000

The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

The instrument ignores any data received from a test or scan for the tier that has sustained the error. Data received from good tiers is retained. All stations in the bad tier are automatically blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Communications problems in instrument.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E18 Normalizer alert – row averages

There is a problem with normalizer panel data. The alert is reported in the following format:

E18	00000001
	00000000
	00000000
	00000000

The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

Normalizer errors are detected during a test cycle. Test data received may be discarded or retained, depending on the error subcode (below). All available stations in the bad tier are automatically blocked; ongoing stations become blocked as testing completes or panels are removed.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- 00000001, 00000002, 00000004, 00000008 – Light source adjustment will be performed when there are no ongoing panels
- 00000010, 00000020, 00000040, 00000080 – All tests are aborted in the affected tier; light source adjustment will be performed when all panels complete or are removed

CORRECTIVE ACTION(S)

- ☐ Automatic light source adjustment will be performed when there are no ongoing panels in the instrument. Do not enter any new panels or move any ongoing panels into any tier reporting this error (i.e., with a subcode other than 00000000). The instrument may be able to correct this error via the light source adjustment process.

E19 Passed read position

One of the tiers passed the position where data scanning must occur. The alert is reported in the following format:

```
E19      00000000
          00000000
          00000100
          00000000
```

The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Optical or mechanical problem with instrument.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E20 Barcode scanner not communicating

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- Instrument cannot communicate properly with a barcode scanner.

00000001 – Tier A scanner

00000002 – Tier B scanner

00000004 – Tier C scanner

00000008 – Tier D scanner

00000010 – External scanner

CORRECTIVE ACTION(S)

- ☐ No testing occurs on the tier experiencing this error, and all the tier stations are automatically blocked. Instrument must be rebooted before the alert can be cleared.

E21 Level 2 rotor step warning

The instrument has sensed that rotor rotation was not ideal but still acceptable. The alert is reported in the following format:

```
E21      00000000
          00000100
          00000000
          00000000
```

The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Likely impediment to or problem with rotor rotation mechanism. Check for a panel protruding from its carrier, an improperly seated closure, or user-applied label peeling off.

CORRECTIVE ACTION(S)

- ☐ If no obvious visible cause for the error exists, contact your local BD representative.

E22 Level 3 rotor step error

The instrument has sensed that rotor rotation was out of specification. The alert is reported in the following format:

```
E22      00000100
          00000000
          00000000
          00000000
```

The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Likely impediment to or problem with rotor rotation mechanism. Check for a panel or panel closure protruding from the carrier or from the panel itself, or for a detached label.

CORRECTIVE ACTION(S)

- ☐ If no obvious visible cause for the error exists, contact your local BD representative.

E23 Normalizer alert – highest/lowest well reading

There is a problem with normalizer panel data. The alert is reported in the following format:

```
E23      00000001
          00000000
          00000000
          00000000
```

The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

Normalizer errors are detected during a test cycle. The instrument ignores any data received from the test for the tier that has sustained the error. Data received from good tiers is retained. All stations in the bad tier are automatically blocked.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- Problem with normalizer panel data.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

E25 UV bulb error

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

Either:

- 00000001 The rear ultraviolet bulb
- 00000002 The front ultraviolet bulb

has burned out and the system cannot switch operation to alternate bulb. The system tried to switch bulbs as a result of an automatic light source adjustment.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative for replacement of the UV bulb(s).

E30 Normalizer Expiration Alert

The instrument issues a weekly alert beginning when Normalizer panel expiration is 60 days away, which progresses to a daily alert when expiration is 30 days away. The alert is reported in the following format:

E30	00000001
	00000001
	00000001
	00000001

The first line of the error represents Normalizers in Tier A, the second line represents Normalizers in Tier B, the third line represents Normalizers in Tier C, and the fourth line represents Normalizers in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- | |
|---|
| 00000001 |
| • Normalizer panel expiration date is between 60 and 30 days away. |
| 00000002 |
| • Normalizer panel expiration date is less than 30 days away (but has not expired). |

CORRECTIVE ACTION(S)

- | |
|---|
| 00000001, 00000002 |
| ☐ Schedule Normalizer panel replacement for the affected tiers before they expire (expiration date is shown on Daily Instrument report). |

E31 Normalizer Expired

Normalizer panel expiration date has passed. Stations in all affected tiers are blocked. Ongoing panel testing is aborted. The alert is reported in the following format:

E31	00000001
	00000001
	00000001
	00000001

The first line of the error represents Normalizers in Tier A, the second line represents Normalizers in Tier B, the third line represents Normalizers in Tier C, and the fourth line represents Normalizers in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- | |
|--|
| 00000001 |
| • Normalizer panel expiration date has passed. |

CORRECTIVE ACTION(S)

- | |
|---|
| 00000001 |
| ☐ Schedule Normalizer panel replacement immediately. |

Update error

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

00000001

- During a “Restore User Data” operation, an error occurred.

00000002

- During an “Install Panel Configuration” operation, an error occurred.

00000004

- During a “Phoenix Update Disk” installation, the number of breakpoints or breakpoint differences exceeded 1200.

00000008

- During a “Phoenix Update Disk” installation, an error occurred.

00000010

- During a “Restore Panel Lot Definition” operation, an error occurred.

CORRECTIVE ACTION(S)

00000001, 00000002, 00000008, 00000010

- ☐ Return to the operation specified above and repeat the process.

00000004

- ☐ Check the breakpoints in the Custom Interpretation Rule Set display. If the number of breakpoints exceeded 1200, breakpoints beyond 1200 were removed. These breakpoints would have appeared at the bottom of the list.

E40 IP Address Network error

The instrument was unable to obtain necessary network information.

ALERT TYPE – 00000001 – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen 00000010 – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- 00000001 – Attempt to connect too many instruments.
- 00000010 – BD EpiCenter system is not running, or hardware problem between systems, such as disconnected cable between instrument and EpiCenter workstation, bad NIC (Network Interface Card) in workstation or instrument, etc.

CORRECTIVE ACTION(S)

- ☐ 00000001 – Do not attempt to connect more than 28 of the same instrument.
- ☐ 00000010 – Start EpiCenter system. Look for obvious source of problem such as disconnected cable. If no obvious source exists, contact your local BD representative.

E44 LIS communications error

The instrument was unable to communicate with the attached LIS system.

ALERT TYPE – Persistent – Cleared when a.) the instrument considers the sub-code reported cleared, and, b.) the alert code has been displayed on the System Alerts screen.

POSSIBLE CAUSE(S)

- 00000001 – LIS system is not running, or hardware problem between systems, such as disconnected cable between instrument and LIS, etc.
- 00000004 – LIS queue is full.

CORRECTIVE ACTION(S)

- ☐ 00000001 – Start LIS system. Look for obvious source of problem such as disconnected cable. If no obvious source exists, contact your local BD representative.
- ☐ 00000004 – Wait until LIS processes some of the pending transmissions from the Phoenix system before attempting to transfer additional ones.

E50 Internal software error

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- System encountered a software error.

CORRECTIVE ACTION(S)

- ☐ Write data to disk (see Section 7.3) and call your local BD representative.

E99 Memory Error

ALERT TYPE – Self-Clearing – cleared when the alert code has been displayed on the System Alerts screen

POSSIBLE CAUSE(S)

- The instrument software encountered a lack of memory.

CORRECTIVE ACTION(S)

- ☐ Contact your local BD representative.

W100 Invalid Panel Sequence Number

POSSIBLE CAUSE(S)

- The panel sequence number typed or scanned does not meet the required number format.

CORRECTIVE ACTION(S)

- ☐ Panel sequence number barcode is located at the top of the reaction side of the panel. Scan or type in the correct panel sequence number. Check panel carton for panel update barcodes; if they are present, refer to Section 4.15 – Installing New Panel Formats, for information on how to proceed.

W101 Missing Panel Sequence Number

POSSIBLE CAUSE(S)

- An operation (find, delete, print, etc.) was attempted and the value in the panel **Sequence #** field was blank.

CORRECTIVE ACTION(S)

- ☐ Type in or scan the correct panel sequence number before attempting the operation again.

W102 Unknown Panel Sequence Number

POSSIBLE CAUSE(S)

- An attempt was made to save or select a panel sequence number that is not in the Phoenix database.

CORRECTIVE ACTION(S)

- ☐ Verify that the correct panel sequence number was entered. If error recurs, the panel may need to be logged in if it is a new panel. Older, completed, finalized panels eventually age out of the database.

W200 Invalid Accession Number

POSSIBLE CAUSE(S)

- An attempt was made to save, find, or print information for an invalid accession number.

CORRECTIVE ACTION(S)

- ☐ Enter a valid accession number, up to 20 characters excluding : * ? [] ! # |

W201 Missing Accession Number

POSSIBLE CAUSE(S)

- An operation (find, save, print, etc.) was attempted and the value in the **Accession #** field was invalid or blank. This can include: a record with a valid panel **Sequence #** and an **Isolate #** greater than 1 with no **Accession #**; an orphan panel with just an **Isolate #**; trying to change a saved record to a blank **Accession #**; etc.

CORRECTIVE ACTION(S)

- ☐ Type in or scan the correct accession number before attempting the operation again.

W202 Unknown Accession Number

POSSIBLE CAUSE(S)

- An attempt was made to find data or print a report for an accession number that is not located in the Phoenix database.

CORRECTIVE ACTION(S)

- ☐ Verify that the correct accession number was entered. If error recurs, the record may need to be logged in if it is new. Older, completed, finalized records eventually age out of the database.

W203 Missing Isolate Number

POSSIBLE CAUSE(S)

- An attempt was made to save a record when an accession number is present without an isolate number.

CORRECTIVE ACTION(S)

- ☐ To save a record, if you enter an accession number, you must also enter an isolate number.

W204 Missing Media Type

POSSIBLE CAUSE(S)

- An attempt was made to save a Yeast ID panel during Panel Login and no Media is specified.

CORRECTIVE ACTION(S)

- ☐ To save a record, if you enter a Yeast ID panel sequence number, you must select a media type.

W300 Missing Test Strain for a QC Panel

POSSIBLE CAUSE(S)

- An attempt was made to save a QC panel without an organism ID (**Test Strain**).

CORRECTIVE ACTION(S)

- ☐ You must select an organism (**Test Strain**) to save a QC panel. Advance to the **Test Strain** field with the "tab" soft key and press DOWN ARROW to drop down a box listing the available test strains. Use UP ARROW or DOWN ARROW to highlight the desired organism. Press the "select" soft key to select the highlighted organism.

W301 Missing Tech ID for a QC Panel

POSSIBLE CAUSE(S)

- An attempt was made to save a QC panel without a Tech ID.

CORRECTIVE ACTION(S)

- ☐ You must enter a Tech ID to save a QC panel. Advance to the **Tech ID** field with the "tab" soft key and enter a Tech ID, up to 3 alphanumeric characters excluding: * ? [] ! # |

W302 Invalid Tech ID

POSSIBLE CAUSE(S)

- An attempt was made to save a QC panel with an invalid Tech ID.

CORRECTIVE ACTION(S)

- ☐ Enter up to 3 alphanumeric characters excluding: * ? [] ! # | for the Tech ID.

W303 Invalid Lot Number**POSSIBLE CAUSE(S)**

- An attempt was made to save a record or print a report for an invalid panel lot number, ID broth lot number, AST broth lot number, or indicator lot number.

CORRECTIVE ACTION(S)

- ☐ Enter the correct lot number. Lot numbers can be up to 7 characters. The lot number is shown on the item carton.

W304 Missing Lot Number**POSSIBLE CAUSE(S)**

- An attempt was made to save or print a report for a QC panel without a panel lot number, or an expiration date has been entered for one of the optional lot number fields (ID Broth, AST Broth, Indicator) with no corresponding lot number.

CORRECTIVE ACTION(S)

- ☐ You must enter a panel lot number to save a QC panel. Advance to the **Panel Lot #** field with the "tab" soft key and enter a lot number, up to 7 digits. The lot number is shown on the panel carton.
- ☐ If you enter an expiration date for one of the optional lot number fields, you must also enter the corresponding lot number.

W305 ID or AST Must be Enabled**POSSIBLE CAUSE(S)**

- An attempt was made to save a panel or QC panel, and both the ID and AST checkboxes are unchecked.

CORRECTIVE ACTION(S)

- ☐ The ID/AST checkboxes are set to checked or unchecked based on the type of panel, according to the panel sequence number. Unchecking either checkbox disables that type of testing. At least one checkbox must be checked for the panel record to be saved and for testing to occur.

W401 Can Not Perform Panel Locate**POSSIBLE CAUSE(S)**

- The panel you are trying to locate is not in the instrument, or you cannot access the instrument using the Normalizer Panel Replacement function (e.g., because there are ongoing panels).

CORRECTIVE ACTION(S)

- ☐ Verify that the correct panel sequence number is entered or that the instrument can be accessed (for Normalizer Panel Replacement activity, ongoing panels render the instrument inaccessible for this activity).

W500 Upgrade Error

POSSIBLE CAUSE(S)

- The floppy disk in the disk drive does not contain a readable, same, or newer version of instrument software. Causes include: older version of instrument software on floppy disk; floppy disk was removed before the update completed; floppy disk contains a corrupted or missing file, etc.

CORRECTIVE ACTION(S)

- ☐ Verify that the correct version (current or later) is printed on the floppy disk label. Attempt the update operation again. If error recurs, contact your local BD representative.

W501 Duplicate Instrument Number

POSSIBLE CAUSE(S)

- An attempt was made to modify the instrument number to be the same as another instrument. (Occurs only when the EpiCenter system is attached.)

CORRECTIVE ACTION(S)

- ☐ Instrument numbers must be unique. Enter a unique instrument number.

W502 Diskette Error

POSSIBLE CAUSE(S)

- The floppy disk drive cannot read a disk. Causes include: no floppy disk in the drive; unformatted disk in the drive; floppy disk has been improperly inserted in the disk drive; floppy disk is or becomes full; floppy disk is write-protected; floppy disk was removed before the operation completed; defective floppy disk drive.

CORRECTIVE ACTION(S)

- ☐ Insert floppy disk with the shutter first and the label facing upward. For save operations, the floppy disk must be formatted, not write-protected, contain sufficient free space for the information being saved, and be left in the drive until the operation completes.

W503 Duplicate Rule

POSSIBLE CAUSE(S)

- In defining custom rules, the data contained in the **Antimicrobial**, **Org(anism) Group**, and **Organism** fields is the same for the rule being saved as a different rule already defined in the Rule Set.

CORRECTIVE ACTION(S)

- ☐ One of the parameters (Antimicrobial, Organism Group, or Organism) must be unique for the record to be saved.

W504 Can Not Change Parameter / Ongoing Panels**POSSIBLE CAUSE(S)**

- An attempt was made to save changes to the **Current Date** or **Current Time** field in Instrument Configuration and Ongoing panels exist in the instrument.

CORRECTIVE ACTION(S)

- ☐ Wait until all panels in the instrument are complete before attempting to perform the operation.

W505 Ineligible For Usage Modification**POSSIBLE CAUSE(S)**

- An attempt was made in Maintenance to Modify Panel Usage and the panel is ineligible for modification. Only combination panels with both ID and AST enabled, that are NOT complete, are eligible for modification.

CORRECTIVE ACTION(S)

- ☐ If incorrect panel Sequence number was entered, enter correct number. Otherwise, panel usage cannot be modified.

W506 Instrument Not Idle**POSSIBLE CAUSE(S)**

- An attempt was made to read light source files or save the syslog from/to floppy disk or to the network and the instrument was not in an idle state (i.e., it was testing, was in the first 30 seconds of an inventory, warming up, etc.), or the carousel was jammed.

CORRECTIVE ACTION(S)

- ☐ Wait until instrument is idle (or correct carousel jam) before performing this activity.

W507 Invalid Password**POSSIBLE CAUSE(S)**

- An invalid (or no) password was entered to access a Configuration display.

CORRECTIVE ACTION(S)

- ☐ Enter the correct password (ECOLI) and press the "perform action" soft key.

W510 Invalid Time Zone**POSSIBLE CAUSE(S)**

- In Instrument Configuration, an attempt was made to save a Time Zone Offset greater than +12 hours or less than -12 hours.

CORRECTIVE ACTION(S)

- ☐ Valid Time Zone Offsets are between -12:00 and +12:00 hours.


W600 Panel Can Not Be Deleted

POSSIBLE CAUSE(S)

- An attempt was made to delete a panel that is still physically in the instrument or a non-pending QC panel that is not physically in the instrument (if QC Lot Support is enabled).

CORRECTIVE ACTION(S)

- ☐ Remove the panel from the instrument before attempting to delete it. You can use the “find panel” soft key (below) to make the carousel present the panel when the door is opened.

- ☐ If you know you have already removed the panel, then close the instrument door and allow the instrument to complete an inventory scan before attempting to delete the panel.

W601 Cannot Finalize

POSSIBLE CAUSE(S)

- An attempt was made to finalize a panel either whose status is not Removable or that has an “unignored” Needs Attention flag set.

CORRECTIVE ACTION(S)

- ☐ If panel is Ongoing, wait until it becomes Removable before attempting to finalize it. If panel has a Needs Attention flag set, resolve the reason before attempting to finalize the panel.

W602 Invalid Date

POSSIBLE CAUSE(S)

- An attempt was made to save a record containing, or to print a report specifying, an invalid date. Or, a panel lot barcode containing an invalid expiration date was scanned.

CORRECTIVE ACTION(S)

- ☐ Refer to **Section 5 – Reference** for information on proper date entry.

W603 Database Underwrite

POSSIBLE CAUSE(S)

- An attempt was made to save displayed information, but the database was modified by another activity, such as an inventory scan, download from the EpiCenter or LIS system, etc.

CORRECTIVE ACTION(S)

- ☐ Recall the information, make the desired changes again, and save the information.

W604 Panel Still Testing**POSSIBLE CAUSE(S)**

- An attempt was made to save a Final ID to an ID or Combination panel that has not reported an Instrument organism ID.
- An attempt was made to perform a Read Light Source from Disk maintenance operation.

CORRECTIVE ACTION(S)

- ☐ A Final ID cannot be selected and saved until the instrument ID results are obtained. Wait until the instrument calculates an ID before attempting to select a Final ID.
- ☐ Wait until all panels in the instrument are complete before attempting to perform the operation.

W605 Invalid Date Range**POSSIBLE CAUSE(S)**

- In Reports, a date greater than the current system date was specified, or a date range exceeding 30 days was entered.

CORRECTIVE ACTION(S)

- ☐ Specify only dates UP TO the current system date.
- ☐ Do not exceed 30 days in date range entries.

W606 Cannot Modify Data for a Finalized Panel**POSSIBLE CAUSE(S)**

- Information for a panel that has been finalized has been recalled and modified. Changes to finalized panels cannot be saved.

CORRECTIVE ACTION(S)

- ☐ If information for a finalized panel must be modified, recall the panel and "un-finalize" it by unchecking the Finalize checkbox. Make the required changes and save the information.

W607 Missing Date**POSSIBLE CAUSE(S)**

- An attempt was made to save a panel without a panel lot number expiration date or without expiration dates for other lot number fields IF panel lot numbers have been entered.

CORRECTIVE ACTION(S)

- ☐ A valid panel lot number expiration date must be entered for panels. If other lot number fields are completed, their corresponding expiration dates must be entered.

W608 Improper Barcode Scan

POSSIBLE CAUSE(S)

- The barcode number scanned is too long for the current active field, or is not the correct type of barcode for the field (e.g., scanning a lot number in a non-lot number field).

CORRECTIVE ACTION(S)

- ☐ Verify which field is currently active and check what type of barcode you are trying to scan. Refer to **Section 5 – Reference** for information on field lengths.

W609 Database Full

POSSIBLE CAUSE(S)

- An attempt was made to log in a new panel and 100 panels exist in the instrument database with a pending record status, with no panels eligible for deletion.

CORRECTIVE ACTION(S)

- ☐ Pending panels should be placed in the instrument and allowed to complete an inventory scan before you attempt to log in any new panels.

W700 Invalid LIS Code

POSSIBLE CAUSE(S)

- An invalid LIS code was entered in the Organism Configuration or Antimicrobial Configuration display.

CORRECTIVE ACTION(S)

- ☐ LIS Codes can be up to 20 alphanumeric characters.

W701 Duplicate LIS Code

POSSIBLE CAUSE(S)

- A duplicate LIS code was entered in the Organism Configuration or Antimicrobial Configuration display.

CORRECTIVE ACTION(S)

- ☐ LIS Codes must be unique in the system.

W702 Invalid LIS Association

POSSIBLE CAUSE(S)

- During Panel Login, when associating downloaded LIS Order records, an attempt was made to associate an Order record containing an organism ID to an ID or Combination panel, or to associate an Order record containing an organism ID to a QC panel, or to associate an Order record containing a Test Strain to a non-QC panel.

CORRECTIVE ACTION(S)

- ☐ You cannot associate an organism ID to a Phoenix ID or Combination panel, nor to a QC panel; and you cannot associate a Test Strain to a non-QC panel.

W800 Panel Lot Undefined

POSSIBLE CAUSE(S)

- An attempt was made to log in a QC panel but the panel sequence number is not within the range of a defined/saved Panel Lot.

CORRECTIVE ACTION(S)

- ☐ When QC Lot Support is enabled in Configuration, all QC panels must belong to a defined/saved panel lot. Clinical panels may be from an undefined panel lot and still be logged in, though such panels generate a Panel Lot Undefined Needs Attention condition.

W801 Panel Lot Range Incomplete

POSSIBLE CAUSE(S)

- During Panel Lot Definition, an attempt was made to save a panel lot but one or both panel sequence number fields were blank.

CORRECTIVE ACTION(S)

- ☐ Both a starting sequence number (lowest) and ending sequence number (highest) must be scanned for any panel lot to be saved. The sequence numbers can be scanned in any order, but both **must** be scanned. You cannot type in a sequence number on the Panel Lot Definition display.

W802 Panel Lot Range Invalid

POSSIBLE CAUSE(S)

- During Panel Lot Definition, an attempt was made to save a panel lot but the two panels scanned were different panel types.

CORRECTIVE ACTION(S)

- ☐ Both panels must be the same type for any panel lot to be saved.

W803 Panel Lot Range Conflict

POSSIBLE CAUSE(S)

- During Panel Lot Definition, an attempt was made to save a panel lot but one or both panel sequence numbers conflicted with an existing panel lot definition.

CORRECTIVE ACTION(S)

- ☐ A panel sequence number can only belong to one defined/saved panel lot. Contact your local BD representative.

7.3 Save System Data to Disk Function

Under certain circumstances, BD may advise you to save (write) system data to a floppy disk. These circumstances include some error conditions and system malfunctions. The Save System Data function is NOT a backup and cannot be restored by the user. To save data to disk, follow the steps below:

- 1 Press the "configuration" soft key.

- 2 The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.

- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE SYSTEM DATA TO DISK."
- 5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.

- 6 Press the "perform action" soft key.

- 7 An hourglass icon appears while the information is being written to the disk.

- 8 If additional disks are required, the system displays the following icon:

Remove the disk and insert another blank, formatted, write-enabled floppy disk in the disk drive.

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button. If there are multiple disks, be sure to number them sequentially and include the date and your facility name on the label.

7.4 Save Syslog to Disk Function

Under certain circumstances, BD may advise you to save (write) the syslog to a floppy disk. These circumstances include some error conditions and system malfunctions. The Save Syslog function copies the system event log, which contains logged system messages about various system, instrument, and communications events. The instrument door must be closed, and the instrument must be idle to save the syslog.

- 1 Press the "configuration" soft key.

- 2 The password window appears. To access configuration functions, enter the appropriate password then press the "perform action" soft key.

- 3 From the main Configuration display, press the "maintenance" soft key.


- 4 Press the UP ARROW or DOWN ARROW to highlight the selection, "SAVE SYSLOG TO DISK."
5 Insert a blank, formatted, write-enabled floppy disk in the disk drive.
6 Press the "perform action" soft key.

- 7 The system asks if you are sure you want to save the syslog. Press the "confirm" soft key to proceed or "cancel" to cancel the operation.

OK

or

- 8 An hourglass icon appears while the information is being written to the disk, as well as the message "Saving Data."

- 9 When the activity is complete, the "activity completed" tone sounds and the hourglass icon disappears. Make sure the disk indicator is off, then remove the floppy disk by pressing the disk eject button.

Performance and Limitations

8.1 Performance Characteristics

Definitions:

Essential Agreement (EA): Essential agreement occurs when the MIC of the Phoenix system and the CLSI Reference Broth Microdilution are identical or within ± 1 dilution of each other.

Category Agreement (CA): Category agreement occurs when the Phoenix system results are within ± 1 dilution from the CLSI Reference Broth Microdilution with respect to the CLSI categorical interpretative criteria.

8.1.1 Gram Negative Performance

Gram Negative Identification

In two internal studies, the performance of the Phoenix Gram Negative Identification was evaluated. The 0.5 inoculum density configuration and the 0.25 inoculum density configuration were tested with 721 strains (0.5) and 784 strains (0.25) respectively. Enteric and non-enteric isolate results were evaluated against commercial and non-commercial methods.

The Phoenix Gram Negative Identification performance is outlined below:

	Inoculum Density (McFarland)	Agreement	No Agreement	No ID
Genus Level	0.5	97.3%	2.1%	0.6%
	0.25	98.5%	1.0%	0.5%
Species Level	0.5	95.6%	3.6%	0.8%
	0.25	98.1%	1.4%	0.5%

Gram Negative Susceptibility

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Amikacin	AN	974	94.3	974	99.2
Amoxicillin	AMX	725	95.9	839	96.9
Amoxicillin-Clavulanate	AMC	636	97.8	749	97.1
Amoxicillin-Clavulanate (f)	AXC	364	96.4	411	97.6
Ampicillin	AM	639	97.5	753	98.0
Ampicillin-Sulbactam	SAM	848	96.8	962	96.9
Ampicillin-Sulbactam (f)	SXA	1682	92.3	1682	94.5
Arbekacin	ARB	2083	94.1	2083	99.2
Aztreonam	ATM	848	96.7	967	97.3
Cefazolin	CZ	634	97.0	752	97.1
Cefdinir	CDR	633	92.7	747	94.5
Cefditoren	CDN	1270	97.2	1270	98.4
Cefepime	FEP	1384	98.0	1384	96.7
Cefetamet	CAT	629	96.0	746	97.7
Cefixime	CFM	1454	96.7	1454	98.1
Cefmetazole	CMZ	608	94.7	718	96.4
Cefoperazone	CFP	854	95.1	972	97.7
Cefoperazone-Sulbactam	SCP	1921	93.1	1921	96.0
Cefotaxime	CTX	849	96.9	970	97.2
Cefotetan	CTT	627	96.2	748	96.4
Cefotiam	CFT	1428	94.9	1428	97.6
Cefoxitin	FOX	628	96.7	748	97.6
Cefpirome	CPO	846	96.2	964	97.0
Cefpodoxime	CPD	616	94.5	737	97.0
Cefsulodin	CFS	55	98.2	56	98.2
Ceftazidime	CAZ	848	95.9	969	97.5
Ceftibuten	CTB	590	90.5	708	95.9
Ceftizoxime	ZOX	854	97.3	971	97.2
Ceftriaxone	CRO	2416	96.1	2416	91.6
Cefuroxime	CXM	623	97.0	744	98.5

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Cephalexin	CN	319	97.8	366	96.5
Cephalothin	CF	613	98.0	731	98.5
Chloramphenicol	C	978	97.7	978	98.7
Ciprofloxacin	CIP	977	99.4	977	99.7
Colistin	CL	467	97.9	467	98.7
Ertapenem*	ETP	1381	98.8	1381	99.6
Fosfomycin	FF	446	93.0	446	96.9
Garenoxacin	GRN	1977	98.1	1977	99.4
Gatifloxacin	GAT	752	99.5	752	99.6
Gemifloxacin	GEM	2096	98.6	2096	99.0
Gentamicin	GM	973	95.8	973	99.8
Imipenem	IPM	975	95.3	975	98.9
Isepamicin	ISP	468	93.8	468	99.4
Kanamycin	K	735	95.8	735	99.9
Levofloxacin	LVX	972	99.6	972	100.0
Lomefloxacin	LOM	976	99.1	976	99.6
Mecillinam	MEC	345	95.7	345	94.8
Meropenem	MEM	969	98.7	969	99.4
Minocycline	MI	2094	93.8	2094	98.3
Moxalactam	MOX	2063	96.9	2063	98.3
Moxifloxacin	MXF	746	98.5	747	99.7
Nalidixic Acid	NA	750	94.0	750	99.1
Netilmicin	NET	974	96.7	974	99.3
Nitrofurantoin	FM	744	98.4	744	98.8
Norfloxacin	NOR	976	98.9	976	99.6
Ofloxacin	OFX	971	99.3	971	99.7
Pefloxacin	PEF	469	98.5	469	99.4
Piperacillin	PIP	860	94.7	973	97.1
Piperacillin-Tazobactam	TZP	856	92.9	970	96.2
Temocillin	TEM	1410	96.1	1410	98.9

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Erythromycin	E	474	95.4	474	95.8
Fosfomycin	FF	472	97.0	472	98.3
Fusidic Acid	FA	477	96.7	477	99.0
Garenoxacin	GRN	1212	97.1	1212	99.0
Gatifloxacin	GAT	481	98.8	481	100.0
Gemifloxacin	GEM	1400	96.2	1400	97.4
Gentamicin	GM	487	98.4	487	99.2
Gentamicin-Syn	GMS			198	97.5
Gentamicin-Syn (SFM)	GMF			87	100.0
Imipenem	IPM	185	100.0	448	97.5
Inducible Macrolide Resistance (iMLSB) Test				295	97.6
Kanamycin	K	463	92.2	464	97.4
Kanamycin Synergy	KS	501	100.0	501	96.6
Levofloxacin	LVX	698	96.9	698	98.6
Lincomycin	L	696	98.3	696	98.7
Linezolid	LZD	474	99.4	474	100
Lomefloxacin	LOM	486	97.7	486	99.2
Meropenem	MEM	185	98.4	447	97.5
Minocycline	MI	1448	98.5	1448	99.2
Moxalactam	MOX	935	95.9	935	96.0
Moxifloxacin	MXF	486	99.0	486	99.6
Mupirocin	MUP	938	96.6	938	99.7
Mupirocin High Level	MUH	968	100.0	968	99.9
Netilmicin	NET	488	93.9	488	96.3
Nitrofurantoin	FM	707	98.3	707	98.3
Norfloxacin	NOR	674	96.0	674	97.8
Ofloxacin	OFX	487	99.2	488	99.8
Oxacillin	OX	449	95.8	449	97.8
Pefloxacin	PEF	706	97.9	706	98.0
Penicillin	P	401	96.0	662	98.0

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Piperacillin	PIP	382	81.7	639	98.0
Piperacillin-Tazobactam	TZP	386	100.0	649	98.3
Pristinamycin	PR	695	93.2	695	95.5
Quinupristin-Dalfopristin	SYN	704	97.2	704	97.6
Rifampin	RA	700	84.0	700	88.3
Streptomycin-Syn	STS			198	99.0
Streptomycin-Syn (SFM)	STF			87	98.8
Teicoplanin	TEC	671	97.6	671	99.4
Telithromycin	TEL	898	97.3	898	98.9
Tetracycline	TE	707	96.3	707	98.6
Ticarcillin	TIC	395	95.4	659	98.2
Ticarcillin-Clavulanate	TIM	396	97.5	656	97.9
Tobramycin	NN	473	97.0	473	97.3
Trimethoprim	TMP	479	95.4	479	97.9
Trimeth-Sulfa (DIN)	STG	645	91.0	646	92.7
Trimethoprim-Sulfamethoxazole	SXT	654	89.5	654	93.0
Vancomycin	VA	702	99.3	702	99.9

* Ability to Detect Resistance Unknown

8.1.3 Streptococci Performance (with the Phoenix SMIC/ID, SMIC Panels)

Streptococci Identification (SMIC/ID only)

In an internal study, the performance of the Phoenix Streptococci Identification was evaluated. Results from 655 isolates were evaluated against commercial and non-commercial methods.

The Phoenix Streptococci Identification performance is outlined below:

	Agreement	No Agreement	No ID
Species Level	96.3%	2.4%	1.2%

Streptococci Susceptibility

The performance of the Phoenix Streptococci AST System was evaluated at four clinical trial sites using clinical isolates. Comparisons were made to AST results generated from reference broth microdilution panels prepared according to CLSI standard guidelines. Discrepant results were arbitrated by duplicate repeat testing in both systems.

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Amoxicillin	AMX	1576	99.2	1576	99.7
Amoxicillin-clavulanate	AMC	1564	97.9	1564	99.3
Ampicillin	AM	1569	98.0	1569	99.4
Cefepime	FEP	1571	97.6	1571	99.7
Cefotaxime	CTX	1578	98.5	1578	99.6
Ceftriaxone	CRO	1579	98.5	1579	99.8
Cefuroxime	CXM	1581	97.7	1581	98.8
Chloramphenicol	C	1587	97.4	1587	99.6
Clindamycin	CC	1591	94.0	1591	97.4
Daptomycin	DAP	1566	93.1	1566	99.9
Ertapenem	ETP	1585	97.9	1585	99.6
Erythromycin	E	1577	93.2	1577	98.2
Garenoxacin*	GRN	1515	95.5	1515	99.7
Gatifloxacin	GAT	1587	95.1	1587	99.5
Gemifloxacin*	GEM	1592	98.9	1592	99.7
Gentamicin-Syn	GMS	1032	99.6	1032	99.6
Imipenem	IPM	1581	98.0	1581	99.6
Kanamycin-Syn	KS	1032	99.6	1032	99.6
Levofloxacin	LVX	1595	97.4	1595	99.8
Linezolid	LZD	1586	97.1	1586	99.4
Meropenem	MEM	1579	97.4	1579	99.7
Moxifloxacin*	MXF	1590	95.8	1590	99.6
Ofloxacin	OFX	1594	98.2	1594	98.7
Penicillin	P	1585	93.9	1585	98.5
Pristinamycin*	PR	1583	94.8	1583	99.8

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Quinupristin-dalfopristin*	SYN	1581	97.3	1581	99.7
Streptomycin-Syn	STS	1031	99.9	1031	99.9
Teicoplanin*	TEC	1593	99.8	1593	99.9
Telithromycin*	TEL	1574	95.9	1574	98.5
Tetracycline	TE	1589	93.3	1589	97.6
Trimethoprim-Sulfamethoxazole	SXT	1008	93.8	1008	98.3
Vancomycin*	VA	1588	98.2	1588	99.9

* Ability to Detect Resistance Unknown

8.1.4 Identification of Yeast Species

The performance of the Phoenix Yeast identification was evaluated across multiple sites using pure colonies isolated from Sabouraud Dextrose Agar (SAB) and Trypticase Soy Agar with 5% Sheep Blood (TSA). Results from 519 (SAB) and 510 (TSA) clinical and challenge isolates were evaluated against conventional and molecular methods.

The Phoenix Yeast identification performance is outlined below:

	Source Media	Agreement	No Agreement	No ID
Genus/Species Level	SAB	95.2%	3.8%	1.0%
	TSA	96.5%	2.7%	0.8%

Additionally, testing was performed at multiple sites to demonstrate reproducibility. The identification results obtained using the Phoenix system were compared with expected results. This performance testing demonstrated inter-site reproducibility of $\geq 95\%$.

8.2 Limitations of the Procedure

See the package insert shipped with the panel for specific organism/antimicrobial limitations.

General

A Gram stain test is required for the selection of the appropriate Phoenix panel types. Accurate identification and/or AST results may not be made without this test.

Use only well-isolated bacterial or yeast colonies from one of the recommended primary isolation media. Use of mixed colonies could result in inaccurate identification and/or AST interpretations.

If the instrument inoculum density (for the panel type being used) is configured to 0.5, an inoculum density of 0.50 – 0.60 McFarland must be met. Only the BBL CrystalSpec or BD PhoenixSpec nephelometer can be used to measure the inoculum density.

If the instrument inoculum density (for the panel type being used) is configured to 0.25, an inoculum density of 0.20 – 0.30 McFarland must be met. Only the BD PhoenixSpec nephelometer can be used to measure the inoculum density for this range.

For identification of yeast, a suspension equivalent of 2.00 - 2.40 McFarland standard must be met and prepared only with the BD PhoenixSpec nephelometer. Use of alternate methods for suspension preparation may cause erroneous identification results.

Phoenix panels can be read only by the Phoenix instrument. Visual interpretation of the Phoenix panels is not possible. Any attempt to manually interpret results from the panel may lead to misidentification and/or inaccurate AST interpretations.

Panels must be placed into the Phoenix instrument within 30 minutes of inoculation.

For the most reliable results, it is recommended that the QC organisms be subcultured at least twice on two consecutive days onto TSA with 5% Sheep Blood before use in the Phoenix system.

Identification

The unique panel environment combined with the shortened incubation time may result in Phoenix panel reactions varying from those obtained using conventional biochemical media.

Antimicrobial Susceptibility Testing

AST Indicator solution may be added to the AST broth tubes and stored in the dark for up to 8 hours prior to use. Use only one free-falling drop of the AST Indicator solution from the dropper bottle. If more than one drop is added inadvertently, discard the tube and use a fresh tube of AST broth.

After the addition of Phoenix AST Indicator Solution to the AST broth tubes, mix by inversion. DO NOT VORTEX. Vortexing may cause air bubbles to form in the AST broth, which can result in inappropriate filling of the Phoenix panel during inoculation.

Because of the low probability of occurrence some organisms included in the ID taxa are not included in the AST database. These organisms will display the message "Organism not included in the AST database, perform alternate method."

For some organism/antimicrobial combinations, the absence of resistant strains precludes defining any result categories other than "susceptible." For strains yielding results suggestive of a "nonsusceptible" category, organism identification and antimicrobial susceptibility test results should be confirmed. Subsequently, the isolates should be saved and submitted to a reference laboratory that will confirm the result using the CLSI reference microbroth dilution method.

The use of CHROMagar Orientation may produce false susceptibility results when testing erythromycin with Gram positive organisms. Antimicrobial susceptibility results should be confirmed using Trypticase Soy Agar with 5% Sheep Blood.

This media type [Chocolate Agar] should *not* be used for Streptococcal identification with SMIC/ID panels. Chocolate Agar may be used for Streptococcal susceptibility testing only.

The use of Columbia Agar with 5% Horse Blood may produce significantly higher MIC for SXT with *Streptococcus* species, which may result in false resistance. Antimicrobial susceptibility test results should be confirmed using Trypticase Soy Agar with 5% Sheep Blood.

NOTES

MIC results for isolates without a genus and species identified in the Phoenix results Final ID field may not be valid.

When changing the Phoenix Final ID from an organism not claimed in the AST Taxa list (see Section F.3) to an organism that is claimed in the AST Taxa, the MIC and Interpretation results will be based on the Final ID that is claimed in the AST Taxa.

Glossary and Abbreviations

absolute threshold

A positivity assessment made by determining if the normalized panel well reading is above (positive) or below (negative) a given fixed threshold.

AST

Antimicrobial Susceptibility Test(ing)

ATCC

American Type Culture Collection

barcode

Machine-readable marking as part of the panel label which uniquely identifies that panel. Contains the product coding and sequence number in Code-128 numeric format.

BP

Blood plate

breakpoint

An interpretation of panel MIC data that which produces Susceptible, Intermediate, or Resistant result classifications. Breakpoints in the MIC data are established by the Clinical and Laboratory Standards Institute (CLSI) and other groups.

caddy

Accessory device used to transport inoculated panels to the instrument for loading.

carousel

The rotating drum which holds the 100 Phoenix panels and positions them for test readings, barcode readings, and loading or unloading.

CLSI

Clinical and Laboratory Standards Institute

CNA

Colistin Nalidixic Acid agar

Col (HRBC)

Columbia agar with Horse Blood

Col (SRBC)

Columbia agar with Sheep Blood

controller

An electronic module that includes a CPU, memory, and interface hardware to control the operation of the instrument.

database

The database is the repository of all information that the instrument collects for the panels that are tested by the instrument. The database is set up to maintain panel results information for the panels tested for 31 days (possibly longer depending on the number of QC panels tested). QC panel results are retained for at least 6 months. The database also maintains the user configurable information for the instrument.

demographic data

Accession information for a panel record

detector

The optical system that measures colorimetric and fluorometric intensities from the reagent reactions in the Phoenix panel wells.

DIN

German Institute for Standardization

DTG

Drug Test Group

EMB

Eosin Methylene Blue agar

end-point results

Results based on well readings obtained after 16 hours of incubation. Microbiological events have reached terminal, or "end-point" values.

error station

Station that has sustained an optical or electro-mechanical error and has been blocked.

EUCAST

European Committee on Antimicrobial Susceptibility Testing

HE

Hektoen Enteric media agar

ID

Microbial Identification

ID/AST combination (combo) panel

The disposable device that contains all reagents needed for both ID and AST.

inoculation station

The inoculation station holds three Phoenix panels at the appropriate angle for optimal fill. The station also holds six broth tubes total, two per organism tested. One tube is for dilution of colony growth for Identification, the other for AST.

instrument test cycle

A complete test of all sample panels located in the carousel, resulting in color and/or fluorescence data values being recorded for each pertinent well of each panel.

MIC

Minimum Inhibitory Concentration; the lowest concentration of an antimicrobial which prohibits continued growth of the tested organism.

normalizer

A reference panel for use in the Phoenix instrument. The Normalizer panel contains a matrix of visible light absorber and fluorescent material in panel-well format, which is used to correct individual well signals for losses occurring in the optical system

orphan

A panel with a valid sequence number, but no associated accession number and isolate number.

panel carrier

The plastic carousel insert, which clips each Phoenix sample panel into place. The carousel contains 104 inserts, 26 in each of four tiers.

panel dataset

Each panel's set of color and fluorescence measurements, the panel's position identifier, test time stamps, and error flags are recorded for each test cycle throughout the test protocol. The test parameters determining each well's results are keyed by the individual barcode label signifying the panel type, and hence the type of test.

PEA

Phenylethyl Alcohol agar

position

The station. The physical location of the Phoenix panel within the instrument. This identifier includes instrument number, carousel tier letter, and numeric position on that tier.

panel presence detection threshold

Each inventory scan or test cycle, the instrument "looks" for both the barcode label of a panel in each carousel position and panel well data, color or fluorescence. If either is detected the instrument declares a panel logically present in that location. A panel without a valid barcode to provide panel type information will not be processed, but will be flagged as a "Needs Attention" candidate.

Phoenix AST indicator

An oxidation-reduction indicator used to signify microbial metabolism in the Phoenix panels. The indicator changes from blue to pink as initial reduction occurs. Further reduction causes the indicator to change from pink to colorless.

"rapid" results

AST Result obtained within 16 hours of panel inoculation.

related panel

Panels with the same accession number and isolate number are related.

Resistance Marker

Condition that is triggered when specific results indicate antimicrobial resistance. The action of some BDxpert rules is to trigger Resistance Markers; other rules may be called as a result of a specific Resistance Marker being triggered.

RGB

Red, Green and Blue. A shorthand representation of the visible light sources / wavelength regions used to interrogate the Phoenix panel.

sample

A specimen contained in a Phoenix panel. In practice, this would be a processed and resuspended dilution of microbiological growth from primary isolation culture in either ID diluent or AST broth which is then poured into the test panel.

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SIR

Susceptible, Intermediate, or Resistant; refers to breakpoint AST categories. See also **Breakpoint**.

standard panel

A special panel containing custom concentrations of various dyes and fluorogenic compounds or materials used to test instrument optical system performance, used internally for engineering, manufacturing and test purposes.

station

The instrument carousel is divided vertically into four tiers (A, B, C, D), each of which holds 26 panels. With one location occupied by a Normalizer panel, 25 locations per tier can accommodate test panels. This means that test panels can populate 100 total locations. Each location is assigned a tier letter and a number to determine the location on the tier. Indicator LEDs located on each panel holder indicate station status (see **Section 3 – Controls and Indicators**).

test sequence number

A count of the number of readings taken by the instrument for a given sample well in a particular test panel, initialized to zero at the time of panel entry.

TSA

Trypticase Soy Agar

XLD

Xylose Lysine Deoxycholate agar

Supplemental ID Test Abbreviations

Abbreviation	Meaning
10C	Growth at 10 degrees Celsius
42C	Growth at 42 degrees Celsius
45C	Growth at 45 degrees Celsius
50C	Growth at 50 degrees Celsius
ACE	MBM + Acetate
ALC	a-Lactose
ANR	Anaerobic growth
ARA	Arabinose
ARG	Arginine
BE	Bile esculin
BSO	Bile solubility
CAMP	CAMP with <i>Staphylococcus aureus</i>
CAT	Catalase
CEL	Cellobiose
CIT	Citrate
COA	Coagulase
DNA	Dnase
ESC	Esculin
FRU	Fructose
GAS	Gas from Glucose
GEL	Gelatin
GLC	Beta glucuronidase
GLU	Glucose
H2S	Hydrogen Sulfide
HGN	Hemolysis - Gram Negative
HGP	Hemolysis - Gram Positive
HIP	Hippurate
IND	Indole
KCN	Growth in KCN
LAC	Lactose
LYS	Lysine
MAC	Growth on MacConkey

Abbreviation	Meaning
MAL	Maltose
MEL	Melibiose
MNS	Mannose
MNT	Mannitol
MOR	Morphology
MOT	Motility
MR	Methyl Red
NIT	Nitrate
NVR	Novobiocin Resistance
OPS	Optochin susceptibility
ORN	Ornithine
OX	Oxidase
PXR	Polymyxin Resistance
PXS	Polymyxin Susceptibility
PYR	Pyrrolidonyl arylamidase
RAF	Raffinose
SBT	Sorbitol
SLT	Growth in 6.5% NaCl
SOR	Sorbose
SUC	Sucrose
TRE	Trehalose
URE	Urea
VAN	Vancomycin
VP	Voges Proskauer
XYL	Xylose
YEL	Yellow/orange pigment

Limited Warranty

This warranty gives you specific legal rights. Additionally, you may have other rights that vary from region to region.

The BD Phoenix™ 100 Automated Microbiology System is warranted to the original purchaser to be free from defects in materials and workmanship for a period of one year following installation. BD's sole responsibility under this warranty shall be to repair or replace any instrument or its components (except for expendable supplies such as printer cartridges, paper, or filters) which under normal operating conditions, prove to be defective within one year of delivery.

BD will furnish new or remanufactured components upon its option. All replacements shall meet new part specifications and shall be warranted as above for the remainder of the one year period. Replaced components become the property of BD.

It is understood that the equipment covered by this Agreement has been installed in accordance with the recommendations and instructions in the *BD Phoenix System User's Manual*.

Any damage to a BD Phoenix system resulting from the insertion or removal of cables that connect this instrument to systems other than those approved or supplied by BD or the failure of the owner to maintain reasonable care and precautions in the operation and maintenance of the system will void this warranty and terminate the obligations of the manufacturer as stated herein.

This warranty is in lieu of all other warranties, whether express or implied, including but not limited to, warranties of merchantability, or fitness for a particular use. In no event will BD be liable for indirect, incidental, special or consequential damages regardless of whether BD has been advised of such.

Replacement Parts

The following items may be ordered by contacting your local BD representative at the numbers in Appendix D.

	Item	Catalog Number
	Barcode Scanner (external)	440503
	Cable, Ethernet	445884
	Cable, Printer	445401
	Cartridge, Laser Printer	445886
	Filter, Air, Electronics Bay, Incubator	447103
	Filter, Air, Source Bay	447101
	Floppy Disks, High Density, Formatted (10)	445847
	Manual, BDPhoenix System User's (ea.)	448021
	Printer, Laser	445882
	Temperature Standard Panel	448984
	Uninterruptible Power Supply (UPS), Europe (1.3 kVA)	445879
	Uninterruptible Power Supply (UPS), Japan (1.3 kVA)	445878
	Uninterruptible Power Supply (UPS), North America (1.3 kVA)	445880

Software Update Log

Whenever you receive a software update, please take a moment to log it below. This can assist you and BD personnel in identifying software revision levels, potential software problems, etc.

Date Received	Software Version	Date Installed	Installed By	Notes

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E.1 Organism Codes, Short Names, Long Names

Organism Long Name	Organism Short Name	Organism Code
<i>Achromobacter denitrificans</i>	Achr. denitrificans	ALCDEN
<i>Achromobacter piechaudii</i>	Achr. piechaudii	ALCPIE
<i>Achromobacter species</i>	Achr. species	ACHRSPE
<i>Achromobacter xylosoxidans ssp xylosoxidans</i>	Achr. xylos. ssp xylosoxidans	ALCXYL
<i>Acinetobacter baumannii</i>	Acinet. baumannii	ACINBAU
<i>Acinetobacter baumannii/calcoaceticus complex</i>	Acinet. baumannii/calco. cplx	ACINBCX
<i>Acinetobacter calcoaceticus</i>	Acinet. calcoaceticus	ACINCAL
<i>Acinetobacter haemolyticus</i>	Acinet. haemolyticus	ACINHAE
<i>Acinetobacter johnsonii</i>	Acinet. johnsonii	ACINJOH
<i>Acinetobacter junii</i>	Acinet. junii	ACINJUN
<i>Acinetobacter lwoffii</i>	Acinet. lwoffii	ACINLWO
<i>Acinetobacter radioresistens</i>	Acinet. radioresistens	ACINRAD
<i>Acinetobacter species</i>	Acinet. species	ACINSPE

Organism Long Name	Organism Short Name	Organism Code
<i>Actinobacillus lignieresii</i>	Actinob. lignieresii	ACTBLIG
<i>Actinobacillus suis</i>	Actinob. suis	ACTBSUI
<i>Actinobacillus ureae</i>	Actinob. ureae	ACTBURE
<i>Aerococcus species</i>	Aeroc. species	AERCSPE
<i>Aerococcus urinae</i>	Aeroc. urinae	AERCURI
<i>Aerococcus viridans</i>	Aeroc. viridans	AERCVIR
<i>Aeromonas allosaccharophila</i>	Aerom. allosaccharophila	AERMALL
<i>Aeromonas caviae</i>	Aerom. caviae	AERMCAV
<i>Aeromonas eucrenophila</i>	Aerom. eucrenophila	AERMEUC
<i>Aeromonas hydrophila</i>	Aerom. hydrophila	AERMHYD
<i>Aeromonas jandaei</i>	Aerom. jandaei	AERMJAN
<i>Aeromonas media</i>	Aerom. media	AERMMED
<i>Aeromonas salmonicida</i>	Aerom. salmonicida	AERMSAL
<i>Aeromonas salmonicida ssp achromogenes</i>	Aerom. salmonic. ssp ach.	AERMSALA
<i>Aeromonas salmonicida ssp masoucida</i>	Aerom. salmonic. ssp mas.	AERMSALM
<i>Aeromonas salmonicida ssp pectinolytica</i>	Aerom. salmonic. ssp pec.	AERMSALPE
<i>Aeromonas salmonicida ssp salmonicida</i>	Aerom. salmonic. ssp sal.	AERMSALSA
<i>Aeromonas salmonicida ssp smithia</i>	Aerom. salmonic. ssp smit.	AERMSALSM
<i>Aeromonas schubertii</i>	Aerom. schubertii	AERMSCH
<i>Aeromonas sobria</i>	Aerom. sobria	AERMVERS
<i>Aeromonas trota</i>	Aerom. trota	AERMTRO
<i>Aeromonas veronii</i>	Aerom. veronii	AERMVERV
<i>Alcaligenes faecalis</i>	Alc. faecalis	ALCFAE
<i>Alloiococcus otitidis</i>	All. otitidis	ALLOTI
<i>Arcanobacterium haemolyticum</i>	Arcan. haemolyticum	ARCAHAE
<i>Arcanobacterium pyogenes</i>	Arcan. pyogenes	ACTMPYO
<i>Aureobasidium pullulans</i>	Aureobs. pullulans	AURPUL
<i>Bacillus cereus</i>	Baci. cereus	BACICER
<i>Bacillus circulans</i>	Baci. circulans	BACICIR

Organism Long Name	Organism Short Name	Organism Code
<i>Bacillus coagulans</i>	Baci. coagulans	BACICOA
<i>Bacillus licheniformis</i>	Baci. licheniformis	BACILIC
<i>Bacillus megaterium</i>	Baci. megaterium	BACIMEG
<i>Bacillus pumilus</i>	Baci. pumilus	BACIPUM
<i>Bacillus subtilis</i>	Baci. subtilis	BACISUB
<i>Bacillus thuringiensis</i>	Baci. thuringiensis	BACITHU
<i>Bergeyella zoohelcum</i>	Ber. zoohelcum	BERZOO
<i>Blastoschizomyces capitatus</i>	Blasto. capitatus	BLACAP
<i>Bordetella bronchiseptica</i>	Bord. bronchiseptica	BORBROS
<i>Brevibacillus brevis</i>	Brevs. brevis	BACIBRE
<i>Brevibacterium species</i>	Brevm. species	BREISPE
<i>Brevundimonas diminuta</i>	Brevu. diminuta	BREUDIM
<i>Brevundimonas vesicularis</i>	Brevu. vesicularis	BREUVES
<i>Burkholderia caryophylli</i>	Burk. caryophylli	BURCAR
<i>Burkholderia cepacia</i>	Burk. cepacia	BURCEP
<i>Burkholderia gladioli</i>	Burk. gladioli	BURGLA
<i>Burkholderia glathei</i>	Burk. glathei	BURGLT
<i>Burkholderia graminis</i>	Burk. graminis	BURGRA
<i>Burkholderia multivorans</i>	Burk. multivorans	BURMUL
<i>Burkholderia phenazinium</i>	Burk. phenazinium	BURPHE
<i>Burkholderia pyrrocinia</i>	Burk. pyrrocinia	BURPYR
<i>Candida albicans</i>	Can. albicans	CANALB
<i>Candida apicola</i>	Can. apicola	CANAPI
<i>Candida boidinii</i>	Can. boidinii	CANBOI
<i>Candida bracarensis</i>	Can. bracarensis	CANBRA
<i>Candida catenulata</i>	Can. catenulata	CANCAT
<i>Candida ciferrii</i>	Can. ciferrii	CANCIF
<i>Candida colliculosa</i>	Can. colliculosa	CANCOL
<i>Candida dubliniensis</i>	Can. dubliniensis	CANDUB

Organism Long Name	Organism Short Name	Organism Code
<i>Candida famata</i>	Can. famata	CANFAM
<i>Candida firmetaria</i>	Can. firmetaria	CANLAM
<i>Candida freyschussii</i>	Can. freyschussii	CANFRE
<i>Candida glabrata</i>	Can. glabrata	TORGLA
<i>Candida guilliermondii</i>	Can. guilliermondii	CANGUI
<i>Candida guilliermondii</i> var <i>membranaefaciens</i>	Can. guillier. var <i>membranaef.</i>	CANGUIM
<i>Candida haemulonii</i>	Can. haemulonii	CANHAE
<i>Candida inconspicua</i>	Can. inconspicua	CANINC
<i>Candida intermedia</i>	Can. intermedia	CANINT
<i>Candida kefyr</i>	Can. kefyr	CANKEF
<i>Candida krusei</i>	Can. krusei	CANKRU
<i>Candida lipolytica</i>	Can. lipolytica	CANLIP
<i>Candida lusitaniae</i>	Can. lusitaniae	CANLUS
<i>Candida magnoliae</i>	Can. magnoliae	CANMAG
<i>Candida melibiosica</i>	Can. melibiosica	CANMEL
<i>Candida membranaefaciens</i>	Can. membranaefaciens	CANMEM
<i>Candida metapsilosis</i>	Can. metapsilosis	CANMET
<i>Candida norvegensis</i>	Can. norvegensis	CANNOR
<i>Candida orthopsilosis</i>	Can. orthopsilosis	CANORT
<i>Candida parapsilosis</i>	Can. parapsilosis	CANPARP
<i>Candida parapsilosis</i> complex	Can. parapsilosis complex	CANPARPX
<i>Candida pararugosa</i>	Can. pararugosa	CANPARR
<i>Candida pelliculosa</i>	Can. pelliculosa	CANPEL
<i>Candida pulcherrima</i>	Can. pulcherrima	CANPUL
<i>Candida rugosa</i>	Can. rugosa	CANRUG
<i>Candida sake</i>	Can. sake	CANSAK
<i>Candida sphaerica</i>	Can. sphaerica	CANSPH
<i>Candida tropicalis</i>	Can. tropicalis	CANTRO
<i>Candida utilis</i>	Can. utilis	CANUTI

Organism Long Name	Organism Short Name	Organism Code
<i>Candida viswanathii</i>	Can. viswanathii	CANVIS
<i>Candida zeylanoides</i>	Can. zeylanoides	CANZEY
<i>Cardiobacterium hominis</i>	Card. hominis	CARHOM
CDC group EF-4a	EF-4a	CDCEF4a
CDC group EF-4b	EF-4b	CDCEF4b
CDC group EO-2	EO-2	CDCEO2
CDC group Vb-3	CDCVb-3	CDCVb3
<i>Cedecea davisae</i>	Ced. davisae	CEDDAV
<i>Cedecea lapagei</i>	Ced. lapagei	CEDLAP
<i>Cedecea neteri</i>	Ced. neteri	CEDNET
<i>Cedecea species 3</i>	Ced. species 3	CEDSPE3
<i>Cedecea species 5</i>	Ced. species 5	CEDSPE5
<i>Cellulomonas turbata</i>	Cell. turbata	OERTUR
<i>Chromobacterium violaceum</i>	Chrom. violaceum	CHROVIO
<i>Chryseobacterium gleum</i>	Chryseob. gleum	CHRBGLE
<i>Chryseobacterium indologenes</i>	Chryseob. indologenes	CHRBIND
<i>Chryseobacterium scophthalmum</i>	Chryseob. scophthalmum	CHRBSCO
<i>Citrobacter amalonaticus</i>	Cit. amalonaticus	CITAMA
<i>Citrobacter braakii</i>	Cit. braakii	CITBRA
<i>Citrobacter farmeri</i>	Cit. farmeri	CITFAR
<i>Citrobacter freundii</i>	Cit. freundii	CITFRE
<i>Citrobacter gillenii</i>	Cit. gillenii	CITSPE10
<i>Citrobacter koseri</i>	Cit. koseri	CITKOS
<i>Citrobacter murlinae</i>	Cit. murlinae	CITSPE11
<i>Citrobacter rodentium</i>	Cit. rodentium	CITSPE9
<i>Citrobacter sedlakii</i>	Cit. sedlakii	CITSED
<i>Citrobacter werkmanii</i>	Cit. werkmanii	CITWER
<i>Citrobacter youngae</i>	Cit. youngae	CITYOU
<i>Comamonas terrigena</i>	Coma. terrigena	COMTER

Organism Long Name	Organism Short Name	Organism Code
<i>Comamonas testosteroni</i>	Coma. testosteroni	COMTES
<i>Corynebacterium amycolatum</i>	Cory. amycolatum	CORAMY
<i>Corynebacterium bovis</i>	Cory. bovis	CORBOV
<i>Corynebacterium diphtheriae</i>	Cory. diphtheriae	CORDIP
<i>Corynebacterium jeikeium</i>	Cory. jeikeium	CORJEI
<i>Corynebacterium kutscheri</i>	Cory. kutscheri	CORKUT
<i>Corynebacterium matruchotii</i>	Cory. matruchotii	CORMAT
<i>Corynebacterium minutissimum</i>	Cory. minutissimum	CORMIN
<i>Corynebacterium propinquum</i>	Cory. propinquum	CORPRO
<i>Corynebacterium pseudodiphtheriticum</i>	Cory. pseudodiphth.	CORPSD
<i>Corynebacterium pseudotuberculosis</i>	Cory. pseudotuberc.	CORPST
<i>Corynebacterium renale</i>	Cory. renale	CORREN
<i>Corynebacterium striatum</i>	Cory. striatum	CORSTR
<i>Corynebacterium ulcerans</i>	Cory. ulcerans	CORULC
<i>Corynebacterium urealyticum</i>	Cory. urealyticum	CORURE
<i>Corynebacterium xerosis</i>	Cory. xerosis	CORXER
<i>Cryptococcus albidus</i>	Cryp. albidus	CRYALB
<i>Cryptococcus gattii</i>	Cryp. gattii	CRYNEOG
<i>Cryptococcus humicolus</i>	Cryp. humicolus	CRYHUM
<i>Cryptococcus laurentii</i>	Cryp. laurentii	CRYLAU
<i>Cryptococcus luteolus</i>	Cryp. luteolus	CRYLUT
<i>Cryptococcus neoformans</i>	Cryp. neoformans	CRYNEO
<i>Cryptococcus terreus</i>	Cryp. terreus	CRYTER
<i>Cryptococcus uniguttulatus</i>	Cryp. uniguttulatus	CRYUNI
<i>Cupriavidus gilardii</i>	Cup. gilardii	RALGIL
<i>Cupriavidus pauculus</i>	Cup. pauculus	CDCIVC2
<i>Delftia acidovorans</i>	Delf. acidovorans	COMACI
<i>Dermabacter hominis</i>	Dermab. hominis	DERBHOM
<i>Dermacoccus nishinomiyaensis</i>	Derm. nishinomiyaen.	MICNIS

Organism Long Name	Organism Short Name	Organism Code
<i>Edwardsiella hoshinae</i>	Ed. hoshinae	EDWHOS
<i>Edwardsiella ictaluri</i>	Ed. ictaluri	EDWICT
<i>Edwardsiella tarda</i>	Ed. tarda	EDWTAR
<i>Eikenella corrodens</i>	Eik. corrodens	EIKCOR
<i>Elizabethkingia meningoseptica</i>	Eliz. meningosept.	CHRBMEN
<i>Empedobacter brevis</i>	Emp. brevis	EMPBRE
<i>Enterobacter aerogenes</i>	Enterob. aerogenes	ENTBAER
<i>Enterobacter amnigenus</i>	Enterob. amnigenus	ENTBAMN
<i>Enterobacter amnigenus</i> biogroup 1	Enterob. amnigen. biogr. 1	ENTBAMN1
<i>Enterobacter amnigenus</i> biogroup 2	Enterob. amnigen. biogr. 2	ENTBAMN2
<i>Enterobacter asburiae</i>	Enterob. asburiae	ENTBASB
<i>Enterobacter cancerogenus</i>	Enterob. cancerogenus	ENTBCAN
<i>Enterobacter cloacae</i>	Enterob. cloacae	ENTBCLO
<i>Enterobacter cloacae</i> ssp <i>dissolvens</i>	Enterob. cloacae ssp dissolven	ENTBDIS
<i>Enterobacter cowanii</i>	Enterob. cowanii	ENTBCOW
<i>Enterobacter gergoviae</i>	Enterob. gergoviae	ENTBGER
<i>Enterobacter hormaechei</i>	Enterob. hormaechei	ENTBHOR
<i>Enterobacter kobei</i>	Enterob. kobei	ENTBKOB
<i>Enterobacter nimipressuralis</i>	Enterob. nimipressuralis	ENTBNIM
<i>Enterobacter sakazakii</i>	Enterob. sakazakii	ENTBSAK
<i>Enterococcus asini</i>	Enteroc. asini	ENTCASI
<i>Enterococcus avium</i>	Enteroc. avium	ENTCAVI
<i>Enterococcus casseliflavus</i>	Enteroc. casseliflavus	ENTCCAS
<i>Enterococcus cecorum</i>	Enteroc. cecorum	ENTCCEC
<i>Enterococcus columbae</i>	Enteroc. columbae	ENTCCOL
<i>Enterococcus dispar</i>	Enteroc. dispar	ENTCDIS
<i>Enterococcus durans</i>	Enteroc. durans	ENTCDUR
<i>Enterococcus faecalis</i>	Enteroc. faecalis	ENTCFAA
<i>Enterococcus faecium</i>	Enteroc. faecium	ENTCFAI

Organism Long Name	Organism Short Name	Organism Code
<i>Enterococcus flavescens</i>	Enteroc. flavescens	ENTCFLA
<i>Enterococcus gallinarum</i>	Enteroc. gallinarum	ENTCGAL
<i>Enterococcus gilvus</i>	Enteroc. gilvus	ENTCGIL
<i>Enterococcus haemoperoxidus</i>	Enteroc. haemoperoxidus	ENTCHAE
<i>Enterococcus hirae</i>	Enteroc. hirae	ENTCHIR
<i>Enterococcus malodoratus</i>	Enteroc. malodoratus	ENTCMAL
<i>Enterococcus moraviensis</i>	Enteroc. moraviensis	ENTCMOR
<i>Enterococcus mundtii</i>	Enteroc. mundtii	ENTCMUN
<i>Enterococcus pallens</i>	Enteroc. pallens	ENTCPAL
<i>Enterococcus pseudoavium</i>	Enteroc. pseudoavium	ENTCPSE
<i>Enterococcus raffinosus</i>	Enteroc. raffinosus	ENTCRAF
<i>Enterococcus ratti</i>	Enteroc. ratti	ENTCRAT
<i>Enterococcus saccharolyticus</i>	Enteroc. saccharolyticus	ENTCSAC
<i>Enterococcus sulfureus</i>	Enteroc. sulfureus	ENTCSUL
<i>Erysipelothrix rhusiopathiae</i>	Ery. rhusiopathiae	ERYRHU
<i>Escherichia blattae</i>	Esch. blattae	ESCBLA
<i>Escherichia coli</i>	Esch. coli	ESCCOL
<i>Escherichia fergusonii</i>	Esch. fergusonii	ESCFER
<i>Escherichia hermannii</i>	Esch. hermannii	ESCHER
<i>Escherichia vulneris</i>	Esch. vulneris	ESCVUL
<i>Ewingella americana</i>	Ew. americana	EWIAME
<i>Exophiala castellani</i>	Exo. castellanii	EXOCAS
<i>Exophiala heteromorpha</i>	Exo. heteromorpha	EXOJEAH
<i>Exophiala jeanselmei</i>	Exo. jeanselmei	EXOJEA
<i>Exophiala lecanii-cornii</i>	Exo. lecanii-cornii	EXOJEAL
<i>Exophiala moniliae</i>	Exo. moniliae	EXOMON
<i>Exophiala species</i>	Exo. species	EXOSPE
<i>Exophiala spinifera</i>	Exo. spinifera	EXOSPI
<i>Gardnerella vaginalis</i>	Gard. vaginalis	GARVAG

Organism Long Name	Organism Short Name	Organism Code
<i>Gemella haemolysans</i>	Gem. haemolysans	GEMHAE
<i>Gemella morbillorum</i>	Gem. morbillorum	GEMMOR
<i>Gemella species</i>	Gem. species	GEMSPE
<i>Geotrichum species</i>	Geo. species	GEOSPE
<i>Globicatella sanguinis</i>	Glob. sanguinis	GLOSAN
<i>Grimontia hollisae</i>	Grim. hollisae	VIBHOL
<i>Hafnia alvei</i>	Haf. alvei	HAFALV
<i>Helcococcus kunzii</i>	Helco. kunzii	HELCKUN
<i>Hortaea werneckii</i>	Hor. werneckii	HORWER
<i>Kingella denitrificans</i>	King. denitrificans	KINDEN
<i>Kingella kingae</i>	King. kingae	KINKIN
<i>Klebsiella granulomatis</i>	Kleb. granulomatis	CALYGRA
<i>Klebsiella oxytoca</i>	Kleb. oxytoca	KLEOXY
<i>Klebsiella pneumoniae ssp ozaenae</i>	Kleb. pneumo. ssp ozae.	KLEPNEO
<i>Klebsiella pneumoniae ssp pneumoniae</i>	Kleb. pneumo. ssp pneu.	KLEPNEP
<i>Klebsiella pneumoniae ssp rhinoscleromatis</i>	Kleb. pneumo. ssp rhino.	KLEPNER
<i>Kloeckera species</i>	Kloeck. species	KLOSPE
<i>Kluyvera ascorbata</i>	Kluyvera ascorbata	KLUASC
<i>Kluyvera cryocrescens</i>	Kluyvera cryocres.	KLUCRY
<i>Kluyvera georgiana</i>	Kluyvera georgiana	KLUGEO
<i>Kluyvera intermedia</i>	Kluyvera intermedia	ENTBINT
<i>Kocuria kristinae</i>	Koc. kristinae	MICKRI
<i>Kocuria rosea</i>	Koc. rosea	MICROS
<i>Kocuria varians</i>	Koc. varians	MICVAR
<i>Kytococcus sedentarius</i>	Kyto. sedentarius	MICSED
<i>Lactococcus garvieae</i>	Lactoc. garvieae	LACCGAR
<i>Lactococcus lactis ssp cremoris</i>	Lactoc. lactis ssp crem.	LACCLACC
<i>Lactococcus lactis ssp hordniae</i>	Lactoc. lactis ssp hord.	LACCLACH
<i>Lactococcus lactis ssp lactis</i>	Lactoc. lactis ssp lactis	LACCLACL

Organism Long Name	Organism Short Name	Organism Code
<i>Lactococcus plantarum</i>	Lactoc. plantarum	LACCPLA
<i>Lactococcus raffinolactis</i>	Lactoc. raffinolactis	LACCRAF
<i>Lactococcus species</i>	Lactoc. species	LACCSPE
<i>Leclercia adecarboxylata</i>	Lec. adecarboxylata	LECADE
<i>Leifsonia aquatica</i>	Leif. aquatica	CORAQU
<i>Leminorella grimontii</i>	Lem. grimontii	LEMGR
<i>Leminorella richardii</i>	Lem. richardii	LEMRI
<i>Leuconostoc citreum</i>	Leu. citreum	LEUCIT
<i>Leuconostoc lactis</i>	Leu. lactis	LEULAC
<i>Leuconostoc mesenteroides ssp cremoris</i>	Leu. mesenter. ssp crem.	LEUMESC
<i>Leuconostoc mesenteroides ssp mesenteroides</i>	Leu. mesenter. ssp mes.	LEUMESM
<i>Leuconostoc pseudomesenteroides</i>	Leu. pseudomesenter.	LEUPSE
<i>Leuconostoc species</i>	Leu. species	LEUSPE
<i>Listeria grayi</i>	Lis. grayi	LISGRA
<i>Listeria innocua</i>	Lis. innocua	LISINN
<i>Listeria ivanovii</i>	Lis. ivanovii	LISIVA
<i>Listeria monocytogenes</i>	Lis. monocytogenes	LISMON
<i>Listeria species</i>	Lis. species	LISSPE
<i>Listeria welshimeri</i>	Lis. welshimeri	LISWEL
<i>Lysinibacillus sphaericus</i>	Lysini. sphaericus	BACISPH
<i>Macrococcus caseolyticus</i>	Macro. caseolyticus	STACAS
<i>Malassezia furfur complex</i>	Mal. furfur complex	MALFURX
<i>Malassezia pachydermatis</i>	Mal. pachydermatis	MALPAC
<i>Malassezia sympodialis</i>	Mal. sympodialis	MALSYM
<i>Mannheimia haemolytica</i>	Mann. haemolytica	PASHAE
<i>Methylobacterium extorquens</i>	Methylob. extorquens	METEXT
<i>Micrococcus luteus</i>	Microc. luteus	MICLUT
<i>Micrococcus lylae</i>	Microc. lylae	MICLYL
<i>Moellerella wisconsensis</i>	Moel. wisconsensis	MOEWIS

Organism Long Name	Organism Short Name	Organism Code
<i>Moraxella (Branhamella) catarrhalis</i>	Morax. (Bran.) cat.	MORABRACAT
<i>Moraxella species</i>	Morax. species	MORASPE
<i>Morganella morganii</i>	Morg. morganii	MORGMOR
<i>Myroides odoratus/odoratimimus</i>	Myr. odoratus/odoratimimus	MYRODAODI
<i>Ochrobactrum anthropi</i>	Och. anthropi	OCHANT
<i>Oerskovia xanthineolytica</i>	Oers. xanthineolytica	OERXAN
<i>Oligella ureolytica</i>	Olig. ureolytica	OLIURO
<i>Oligella urethralis</i>	Olig. urethralis	OLIURT
<i>Paenibacillus alvei</i>	Paen. alvei	PAEALV
<i>Paenibacillus macerans</i>	Paen. macerans	PAEMAC
<i>Pantoea agglomerans</i>	Pan. agglomerans	PANAGG
<i>Pantoea ananatis</i>	Pan. ananatis	ERWANA
<i>Pantoea dispersa</i>	Pan. dispersa	PANDIS
<i>Pantoea stewartii ssp indologenes</i>	Pan. ste. ssp indologenes	PANSTEI
<i>Pantoea stewartii ssp stewartii</i>	Pan. ste. ssp stewartii	PANSTES
<i>Pasteurella aerogenes</i>	Past. aerogenes	PASAER
<i>Pasteurella multocida</i>	Past. multocida	PASMUL
<i>Pasteurella pneumotropica</i>	Past. pneumotropica	PASPNE
<i>Pediococcus acidilactici</i>	Ped. acidilactici	PEDACI
<i>Pediococcus damnosus</i>	Ped. damnosus	PEDDAM
<i>Pediococcus dextrinicus</i>	Ped. dextrinicus	PEDDEX
<i>Pediococcus parvulus</i>	Ped. parvulus	PEDPAR
<i>Pediococcus pentosaceus</i>	Ped. pentosaceus	PEDPEN
<i>Pediococcus species</i>	Ped. species	PEDSPE
<i>Photobacterium damsela</i>	Photob. damsela	PHOBDAM
<i>Pichia angusta</i>	Pic. angusta	PICANG
<i>Pichia burtonii</i>	Pic. burtonii	PICBUR
<i>Pichia farinosa</i>	Pic. farinosa	PICFAR
<i>Plesiomonas shigelloides</i>	Ples. shigelloides	PLESHI

Organism Long Name	Organism Short Name	Organism Code
<i>Pragia fontium</i>	Prag. fontium	PRAFON
<i>Proteus hauseri</i>	Prot. hauseri	PROTHAU
<i>Proteus mirabilis</i>	Prot. mirabilis	PROTMIR
<i>Proteus myxofaciens</i>	Prot. myxofaciens	PROTMYX
<i>Proteus penneri</i>	Prot. penneri	PROTPEN
<i>Proteus vulgaris</i>	Prot. vulgaris	PROTVUL
<i>Prototheca species</i>	Protot. species	PROHSPE
<i>Prototheca wickerhamii</i>	Protot. wick.	PROHWIC
<i>Prototheca zopfii</i>	Protot. zopfii	PROHZOP
<i>Providencia alcalifaciens</i>	Prov. alcalifaciens	PROVALC
<i>Providencia heimbachae</i>	Prov. heimbachae	PROVHEI
<i>Providencia rettgeri</i>	Prov. rettgeri	PROVRET
<i>Providencia rustigianii</i>	Prov. rustigianii	PROVRUS
<i>Providencia stuartii</i>	Prov. stuartii	PROVSTU
<i>Pseudomonas aeruginosa</i>	Pseud. aeruginosa	PSEAER
<i>Pseudomonas alcaligenes</i>	Pseud. alcaligenes	PSEALC
<i>Pseudomonas fluorescens</i>	Pseud. fluorescens	PSEFLU
<i>Pseudomonas luteola</i>	Pseud. luteola	CHRM LUT
<i>Pseudomonas mendocina</i>	Pseud. mendocina	PSEMEN
<i>Pseudomonas monteilii</i>	Pseud. monteilii	PSEMON
<i>Pseudomonas oryzihabitans</i>	Pseud. oryzihabitans	FLAIORY
<i>Pseudomonas pertucinogena</i>	Pseud. pertucinogena	PSEPER
<i>Pseudomonas pseudoalcaligenes</i>	Pseud. pseudoalcaligenes	PSEPSE
<i>Pseudomonas putida</i>	Pseud. putida	PSEPUT
<i>Pseudomonas species</i>	Pseud. species	PSESPE
<i>Pseudomonas stutzeri</i>	Pseud. stutzeri	PSESTU
<i>Pseudomonas veronii</i>	Pseud. veronii	PSEVER
<i>Rahnella aquatilis</i>	Rah. aquatilis	RAHAQU
<i>Ralstonia pickettii</i>	Ral. pickettii	BURPIC

Organism Long Name	Organism Short Name	Organism Code
<i>Ralstonia solanacearum</i>	Ral. solanacearum	BURSOL
<i>Ralstonia species</i>	Ral. species	RALSPE
<i>Raoultella ornithinolytica</i>	Rao. ornithinolytica	KLEORN
<i>Raoultella planticola</i>	Rao. planticola	KLEPLA
<i>Raoultella terrigena</i>	Rao. terrigena	KLETER
<i>Rhizobium radiobacter</i>	Rhizob. radiobacter	AGRRAD
<i>Rhodococcus equi</i>	Rhodoc. equi	RHOCEQU
<i>Rhodotorula glutinis</i>	Rhodot. glutinis	RHOTGLU
<i>Rhodotorula minuta</i>	Rhodot. minuta	RHOTMIN
<i>Rhodotorula mucilaginosa</i> var <i>mucilaginosa</i>	Rhodot. mucilag. var mucilag.	RHOTMUCM
<i>Rothia dentocariosa</i>	Roth. dentocariosa	ROTDEN
<i>Rothia mucilaginosa</i>	Roth. mucilaginosa	STOMUC
<i>Saccharomyces cerevisiae</i>	Sac. cerevisiae	SACCER
<i>Salmonella aberdeen</i>	Salm. aberdeen	SALABE
<i>Salmonella abortus-equi</i>	Salm. abortus-equi	SALABOE
<i>Salmonella adelaide</i>	Salm. adelaide	SALADL
<i>Salmonella aderike</i>	Salm. aderike	SALADR
<i>Salmonella agona</i>	Salm. agona	SALAGO
<i>Salmonella alachua</i>	Salm. alachua	SALALA
<i>Salmonella anatum</i>	Salm. anatum	SALANA
<i>Salmonella arizonae</i>	Salm. arizonae	SALARI
<i>Salmonella avana</i>	Salm. avana	SALAVA
<i>Salmonella bahrenfeld</i>	Salm. bahrenfeld	SALBAH
<i>Salmonella blockley</i>	Salm. blockley	SALBLO
<i>Salmonella bongori</i>	Salm. bongori	SALBON
<i>Salmonella braenderup</i>	Salm. braenderup	SALBRA
<i>Salmonella bredeney</i>	Salm. bredeney	SALBRE
<i>Salmonella bunn</i>	Salm. bunn	SALBUN
<i>Salmonella californica</i>	Salm. californica	SALCAL

Organism Long Name	Organism Short Name	Organism Code
<i>Salmonella carrau</i>	Salm. carrau	SALCAR
<i>Salmonella cerro</i>	Salm. cerro	SALCER
<i>Salmonella champaign</i>	Salm. champaign	SALCHA
<i>Salmonella chittagong</i>	Salm. chittagong	SALCHI
<i>Salmonella choleraesuis</i>	Salm. choleraesuis	SALCHO
<i>Salmonella choleraesuis ssp arizonae</i>	Salm. choleraesuis ssp ariz.	SALCHOA
<i>Salmonella choleraesuis ssp choleraesuis</i>	Salm. choleraesuis ssp chol.	SALCHOC
<i>Salmonella choleraesuis ssp diarizonae</i>	Salm. choleraesuis ssp diariz.	SALCHOD
<i>Salmonella choleraesuis ssp houtenae</i>	Salm. choleraesuis ssp hout.	SALCHOH
<i>Salmonella choleraesuis ssp indica</i>	Salm. choleraesuis ssp indica	SALCHOI
<i>Salmonella choleraesuis ssp salamae</i>	Salm. choleraesuis ssp salamae	SALCHOS
<i>Salmonella cubana</i>	Salm. cubana	SALCUB
<i>Salmonella dakar</i>	Salm. dakar	SALDAK
<i>Salmonella daressalaam</i>	Salm. daressalaam	SALDAR
<i>Salmonella derby</i>	Salm. derby	SALDER
<i>Salmonella dessau</i>	Salm. dessau	SALDES
<i>Salmonella DT</i>	Salm. DT	SALDT
<i>Salmonella dublin</i>	Salm. dublin	SALDUB
<i>Salmonella duesseldorf</i>	Salm. duesseldorf	SALDUE
<i>Salmonella enteritidis</i>	Salm. enteritidis	SALENT
<i>Salmonella fresno</i>	Salm. fresno	SALFRE
<i>Salmonella gallinarum</i>	Salm. gallinarum	SALGAL
<i>Salmonella gaminara</i>	Salm. gaminara	SALGAM
<i>Salmonella give</i>	Salm. give	SALGIV
<i>Salmonella haardt</i>	Salm. haardt	SALHAA
<i>Salmonella hadar</i>	Salm. hadar	SALHAD
<i>Salmonella hamburg</i>	Salm. hamburg	SALHAM
<i>Salmonella hartford</i>	Salm. hartford	SALHAR
<i>Salmonella heidelberg</i>	Salm. heidelberg	SALHEI

Organism Long Name	Organism Short Name	Organism Code
<i>Salmonella illinois</i>	Salm. illinois	SALILL
<i>Salmonella infantis</i>	Salm. infantis	SALINF
<i>Salmonella invernness</i>	Salm. invernness	SALINV
<i>Salmonella java</i>	Salm. java	SALJAA
<i>Salmonella javiana</i>	Salm. javiana	SALJAI
<i>Salmonella kentucky</i>	Salm. kentucky	SALKEN
<i>Salmonella kirkee</i>	Salm. kirkee	SALKIR
<i>Salmonella kunduchi</i>	Salm. kunduchi	SALKUN
<i>Salmonella kvittingfoss</i>	Salm. kvittingfoss	SALKVI
<i>Salmonella lansing</i>	Salm. lansing	SALLAN
<i>Salmonella litchfield</i>	Salm. litchfield	SALLIT
<i>Salmonella liverpool</i>	Salm. liverpool	SALLIV
<i>Salmonella london</i>	Salm. london	SALLON
<i>Salmonella luciana</i>	Salm. luciana	SALLUC
<i>Salmonella manhattan</i>	Salm. manhattan	SALMAN
<i>Salmonella mbandaka</i>	Salm. mbandaka	SALMBA
<i>Salmonella meleagridis</i>	Salm. meleagridis	SALMEL
<i>Salmonella memphis</i>	Salm. memphis	SALMEM
<i>Salmonella michigan</i>	Salm. michigan	SALMIC
<i>Salmonella minneapolis</i>	Salm. minneapolis	SALMIA
<i>Salmonella minnesota</i>	Salm. minnesota	SALMIS
<i>Salmonella montevideo</i>	Salm. montevideo	SALMON
<i>Salmonella muenchen</i>	Salm. muenchen	SALMUE
<i>Salmonella muenster</i>	Salm. muenster	SALMUS
<i>Salmonella newington</i>	Salm. newington	SALNEI
<i>Salmonella newport</i>	Salm. newport	SALNEP
<i>Salmonella nottingham</i>	Salm. nottingham	SALNOT
<i>Salmonella ohio</i>	Salm. ohio	SALOH1
<i>Salmonella onderstepoort</i>	Salm. onderstepoort	SALOND

Organism Long Name	Organism Short Name	Organism Code
<i>Salmonella oranienburg</i>	Salm. oranienburg	SALORA
<i>Salmonella panama</i>	Salm. panama	SALPAN
<i>Salmonella paratyphi A</i>	Salm. paratyphi A	SALPARA
<i>Salmonella paratyphi B</i>	Salm. paratyphi B	SALPARB
<i>Salmonella poona</i>	Salm. poona	SALPOO
<i>Salmonella pullorum</i>	Salm. pullorum	SALPUL
<i>Salmonella quinhon</i>	Salm. quinhon	SALQUI
<i>Salmonella rubislaw</i>	Salm. rubislaw	SALRUB
<i>Salmonella saintpaul</i>	Salm. saintpaul	SALSAI
<i>Salmonella schwarzengrund</i>	Salm. schwarzengrund	SALSCH
<i>Salmonella senftenberg</i>	Salm. senftenberg	SALSEN
<i>Salmonella species</i>	Salm. species	SALSPE
<i>Salmonella tallahassee</i>	Salm. tallahassee	SALTAL
<i>Salmonella thompson</i>	Salm. thompson	SALTHO
<i>Salmonella typhi</i>	Salm. typhi	SALTYP
<i>Salmonella typhimurium</i>	Salm. typhimurium	SALTYPM
<i>Salmonella virginia</i>	Salm. virginia	SALVIR
<i>Salmonella westerstede</i>	Salm. westerstede	SALWES
<i>Salmonella worthington</i>	Salm. worthington	SALWOR
<i>Serratia entomophila</i>	Ser. entomophila	SERENT
<i>Serratia ficaria</i>	Ser. ficaria	SERFIC
<i>Serratia fonticola</i>	Ser. fonticola	SERFON
<i>Serratia grimesii</i>	Ser. grimesii	SERGRI
<i>Serratia liquefaciens</i>	Ser. liquefaciens	SERLIQ
<i>Serratia marcescens</i>	Ser. marcescens	SERMAR
<i>Serratia odorifera</i>	Ser. odorifera	SERODO
<i>Serratia odorifera 1</i>	Ser. odorifera 1	SERODO1
<i>Serratia odorifera 2</i>	Ser. odorifera 2	SERODO2
<i>Serratia plymuthica</i>	Ser. plymuthica	SERPLY

Organism Long Name	Organism Short Name	Organism Code
<i>Serratia proteamaculans ssp proteamaculans</i>	Ser. proteamac. ssp proteam.	SERPROP
<i>Serratia proteamaculans ssp quinovora</i>	Ser. proteamac. ssp quino.	SERPROQ
<i>Serratia rubidaea</i>	Ser. rubidaea	SERRUB
<i>Shewanella algae</i>	Shew. algae	SHEALG
<i>Shewanella putrefaciens</i>	Shew. putrefaciens	SHEPUT
<i>Shigella boydii</i>	Shig. boydii	SHIBOY
<i>Shigella dysenteriae</i>	Shig. dysenteriae	SHIDYS
<i>Shigella flexneri</i>	Shig. flexneri	SHIFLE
<i>Shigella sonnei</i>	Shig. sonnei	SHISON
<i>Shigella species</i>	Shig. species	SHISPE
<i>Sphingobacterium multivorum</i>	Sphingob. multivorum	SPHBMUL
<i>Sphingobacterium spiritivorum</i>	Sphingob. spiritivorum	SPHBSPI
<i>Sphingobacterium thalpophilum</i>	Sphingob. thalpophilum	SPHBTHA
<i>Sphingomonas paucimobilis</i>	Sphingom. paucimobilis	SPHMPAU
<i>Sporobolomyces salmonicolor</i>	Sporobol. salmonicolor	SPOBSAL
<i>Staphylococcus arlettae</i>	Staph. arlettae	STAARL
<i>Staphylococcus aureus</i>	Staph. aureus	STAAUE
<i>Staphylococcus auricularis</i>	Staph. auricularis	STAAUI
<i>Staphylococcus capitis</i>	Staph. capitis	STACAI
<i>Staphylococcus capitis ssp capitis</i>	Staph. capitis ssp capitis	STACAIC
<i>Staphylococcus capitis ssp ureolyticus</i>	Staph. capitis ssp ureolyt.	STACAIU
<i>Staphylococcus caprae</i>	Staph. caprae	STACAP
<i>Staphylococcus carnosus</i>	Staph. carnosus	STACAR
<i>Staphylococcus chromogenes</i>	Staph. chromogenes	STACHR
<i>Staphylococcus cohnii</i>	Staph. cohnii	STACOH
<i>Staphylococcus cohnii ssp cohnii</i>	Staph. cohnii ssp cohnii	STACOHC
<i>Staphylococcus cohnii ssp urealyticum</i>	Staph. cohnii ssp urealyt.	STACOHU
<i>Staphylococcus condimentii</i>	Staph. condimentii	STACON
<i>Staphylococcus delphini</i>	Staph. delphini	STADEL

Organism Long Name	Organism Short Name	Organism Code
<i>Staphylococcus epidermidis</i>	Staph. epidermidis	STAEPI
<i>Staphylococcus equorum</i>	Staph. equorum	STAEQU
<i>Staphylococcus felis</i>	Staph. felis	STAFEL
<i>Staphylococcus fleurettii</i>	Staph. fleurettii	STAFLE
<i>Staphylococcus gallinarum</i>	Staph. gallinarum	STAGAL
<i>Staphylococcus haemolyticus</i>	Staph. haemolyticus	STAHAE
<i>Staphylococcus hominis</i>	Staph. hominis	STAHOM
<i>Staphylococcus hyicus</i>	Staph. hyicus	STAHYI
<i>Staphylococcus intermedius</i>	Staph. intermedius	STAINI
<i>Staphylococcus kloosii</i>	Staph. kloosii	STAKLO
<i>Staphylococcus lentus</i>	Staph. lentus	STALEN
<i>Staphylococcus lugdunensis</i>	Staph. lugdunensis	STALUG
<i>Staphylococcus lutrae</i>	Staph. lutrae	STALUT
<i>Staphylococcus muscae</i>	Staph. muscae	STAMUS
<i>Staphylococcus pasteurii</i>	Staph. pasteurii	STAPAS
<i>Staphylococcus piscifermentans</i>	Staph. piscifermentans	STAPIS
<i>Staphylococcus pulvereri</i>	Staph. pulvereri	STAPUL
<i>Staphylococcus saccharolyticus</i>	Staph. saccharolyticus	STASAC
<i>Staphylococcus saprophyticus</i>	Staph. saprophyticus	STASAP
<i>Staphylococcus schleiferi</i>	Staph. schleiferi	STASCH
<i>Staphylococcus schleiferi ssp coagulans</i>	Staph. schleiferi ssp coagul.	STASCHC
<i>Staphylococcus schleiferi ssp schleiferi</i>	Staph. schleiferi ssp schleif.	STASCHS
<i>Staphylococcus sciuri</i>	Staph. sciuri	STASCI
<i>Staphylococcus simulans</i>	Staph. simulans	STASIM
<i>Staphylococcus succinus ssp casei</i>	Staph. suc. ssp casei	STASUCCA
<i>Staphylococcus succinus ssp succinus</i>	Staph. suc. ssp succinus	STASUCSU
<i>Staphylococcus vitulinus</i>	Staph. vitulinus	STAVIT
<i>Staphylococcus warneri</i>	Staph. warneri	STAWAR
<i>Staphylococcus xylosus</i>	Staph. xylosus	STAXYL

Organism Long Name	Organism Short Name	Organism Code
<i>Stenotrophomonas africana</i>	Sten. africana	STEAFR
<i>Stenotrophomonas maltophilia</i>	Sten. maltophilia	STEMAL
<i>Streptococcus acidominimus</i>	Strep. acidominimus	STRACI
<i>Streptococcus agalactiae</i> (Strep. group B)	Strep. agalactiae (Str. gr. B)	STRAGA
<i>Streptococcus alactolyticus</i>	Strep. alactolyticus	STRALA
<i>Streptococcus anginosus</i>	Strep. anginosus	STRANG
<i>Streptococcus anginosus</i> (previously <i>milleri</i>) group	Strep. anginosus (<i>milleri</i>) gr.	STRANGGR
<i>Streptococcus bovis</i> (Strep. group D)	Strep. bovis (Str. gr. D)	STRBOV
<i>Streptococcus bovis</i> I (Strep. group D)	Strep. bovis I (Str. gr. D)	STRBOVI
<i>Streptococcus bovis</i> II (Strep. group D)	Strep. bovis II (Str. gr. D)	STRBOVII
<i>Streptococcus canis</i>	Strep. canis	STRCAN
<i>Streptococcus constellatus</i>	Strep. constellatus	STRCON
<i>Streptococcus criceti</i>	Strep. criceti	STRCRC
<i>Streptococcus cristatus</i>	Strep. cristatus	STRCRS
<i>Streptococcus downei</i>	Strep. downei	STRDOW
<i>Streptococcus dysgalactiae</i>	Strep. dysgalactiae	STRDYS
<i>Streptococcus dysgalactiae</i> ssp <i>dysgalactiae</i>	Strep. dysgal. ssp dysgal.	STRDYSDY
<i>Streptococcus dysgalactiae</i> ssp <i>equisimilis</i>	Strep. dysgal. ssp equis.	STRDYSEM
<i>Streptococcus dysgalactiae/canis</i>	Strep. dysgal./canis	STRDYSCAN
<i>Streptococcus equi</i>	Strep. equi	STREQU
<i>Streptococcus equi</i> ssp <i>equi</i>	Strep. equi ssp equi	STREQUE
<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>	Strep. equi ssp zooepid.	STREQUZ
<i>Streptococcus equinus</i>	Strep. equinus	STREQN
<i>Streptococcus ferus</i>	Strep. ferus	STRFER
<i>Streptococcus gallolyticus</i> ssp <i>macedonicus</i>	Strep. gallolyti. ssp macedon.	STRMAC
<i>Streptococcus gordonii</i>	Strep. gordonii	STRGOR
<i>Streptococcus</i> group A (<i>Strep. pyogenes</i>)	Strep. group A (Str. pyogenes)	STRGRA
<i>Streptococcus</i> group B (<i>Strep. agalactiae</i>)	Strep. group B (Str. agalact.)	STRGRB

Organism Long Name	Organism Short Name	Organism Code
<i>Streptococcus</i> group C (large colony)	Strep. group C (lg col)	STRGRC
<i>Streptococcus</i> group C/G (large colony)	Strep. group C/G (lg col)	STRGRCG
<i>Streptococcus</i> group CFG (small colony)	Strep. group CFG (sm col)	STRGRCFG
<i>Streptococcus</i> group D (non-enterococcus)	Strep. group D (non-enteroc.)	STRGRDNE
<i>Streptococcus</i> group E	Strep. group E	STRGRE
<i>Streptococcus</i> group F	Strep. group F	STRGRF
<i>Streptococcus</i> group G (large colony)	Strep. group G (lg col)	STRGRG
<i>Streptococcus</i> group L	Strep. group L	STRGRL
<i>Streptococcus hyointestinalis</i>	Strep. hyointestinalis	STRHYO
<i>Streptococcus infantarius ssp coli</i>	Strep. infa ssp coli	STRINACO
<i>Streptococcus infantarius ssp infantarius</i>	Strep. infa ssp infantarius	STRINAIN
<i>Streptococcus infantis</i>	Strep. infantis	STRINF
<i>Streptococcus iniae</i>	Strep. iniae	STRINI
<i>Streptococcus intermedius</i>	Strep. intermedius	STRINR
<i>Streptococcus milleri</i> group	Strep. milleri gr.	STRMILGR
<i>Streptococcus mitis</i>	Strep. mitis	STRMIT
<i>Streptococcus mitis</i> group	Strep. mitis gr.	STRMITGR
<i>Streptococcus mutans</i>	Strep. mutans	STRMUT
<i>Streptococcus mutans</i> group	Strep. mutans gr.	STRMUTGR
<i>Streptococcus oralis</i>	Strep. oralis	STRORA
<i>Streptococcus parasanguinis</i>	Strep. parasanguinis	STRPAR
<i>Streptococcus peroris</i>	Strep. peroris	STRPER
<i>Streptococcus pneumoniae</i>	Strep. pneumoniae	STRPNE
<i>Streptococcus porcinus</i>	Strep. porcinus	STRPOR
<i>Streptococcus pyogenes</i> (Strep. group A)	Strep. pyogenes (Str. gr. A)	STRPYO
<i>Streptococcus ratti</i>	Strep. ratti	STRRAT
<i>Streptococcus salivarius</i>	Strep. salivarius	STRSAL
<i>Streptococcus salivarius</i> group	Strep. salivarius gr.	STRSALGR
<i>Streptococcus sanguinis</i>	Strep. sanguinis	STRSAN

Organism Long Name	Organism Short Name	Organism Code
<i>Streptococcus sanguinis</i> group	Strep. sanguinis gr.	STRSANGR
<i>Streptococcus sobrinus</i>	Strep. sobrinus	STRSOB
<i>Streptococcus suis</i>	Strep. suis	STRSUI
<i>Streptococcus thermophilus</i>	Strep. thermophilus	STRTHE
<i>Streptococcus uberis</i>	Strep. uberis	STRUBE
<i>Streptococcus vestibularis</i>	Strep. vestibularis	STRVES
<i>Streptococcus viridans</i> group	Strep. viridans group	STRVIRGR
<i>Suttonella indologenes</i>	Sutto. indologenes	SUTIND
<i>Tatumella pyseos</i>	Tat. pyseos	TATPTY
<i>Trichosporon asahii</i>	Tric. asahii	TRIASA
<i>Trichosporon asteroides</i>	Tric. asteroides	TRIAST
<i>Trichosporon cutaneum</i>	Tric. cutaneum	TRICUT
<i>Trichosporon inkin</i>	Tric. inkin	TRIINK
<i>Trichosporon loubieri</i>	Tric. loubieri	TRILOU
<i>Trichosporon mucoides</i>	Tric. mucoides	TRIMUC
<i>Trichosporon ovoides</i>	Tric. ovoides	TRIOVO
<i>Trichosporon pullulans</i>	Tric. pullulans	TRIPUL
<i>Vibrio alginolyticus</i>	Vib. alginolyticus	VIBALG
<i>Vibrio cholerae</i>	Vib. cholerae	VIBCHO
<i>Vibrio fluvialis</i>	Vib. fluvialis	VIBFLU
<i>Vibrio metschnikovii</i>	Vib. metschnikovii	VIBMET
<i>Vibrio mimicus</i>	Vib. mimicus	VIBMIM
<i>Vibrio parahaemolyticus</i>	Vib. parahaemolyticus	VIBPAR
<i>Vibrio vulnificus</i>	Vib. vulnificus	VIBVUL
<i>Wangiella dermatitidis</i>	Wan. dermatitidis	WANDER
<i>Weeksella virosa</i>	Week. virosa	WEEVIR
<i>Yersinia aldovae</i>	Yer. aldovae	YERALD
<i>Yersinia bercovieri</i>	Yer. bercovieri	YERBER
<i>Yersinia enterocolitica</i>	Yer. enterocolitica	YERENT

Organism Long Name	Organism Short Name	Organism Code
<i>Yersinia frederiksenii</i>	Yer. frederiksenii	YERFRE
<i>Yersinia intermedia</i>	Yer. intermedia	YERINT
<i>Yersinia kristensenii</i>	Yer. kristensenii	YERKRI
<i>Yersinia mollaretii</i>	Yer. mollaretii	YERMOL
<i>Yersinia pseudotuberculosis</i>	Yer. pseudotuberculosis	YERPSE
<i>Yersinia rohdei</i>	Yer. rohdei	YERROH
<i>Yersinia ruckeri</i>	Yer. ruckeri	YERRUC
<i>Yokenella regensburgei</i>	Yok. regensburgei	YOKREG
<i>Zygosaccharomyces bailli</i>	Zyg. bailli	ZYGBAI

Panel Information

F.1 List of Antimicrobial Agents in Phoenix Panels

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Class	Drug Name ↓	Drug Code	Range	Range	Range
Aminoglycoside	Amikacin	AN	0.5-64	0.5-64	N/A
Beta-lactam	Amoxicillin	AMX	0.5-32	0.25-32	0.125-32
Beta-lactam	Amoxicillin-Clavulanate (f)	AXC	0.5/2-32/2	0.25/2-32/2	N/A
Beta-lactam	Amoxicillin-Clavulanate	AMC	0.5/0.25-32/16	0.25/0.12-32/16	0.125/0.06-32/16
Beta-lactam	Ampicillin	AM	0.5-32	0.06-32	0.06-32
Beta-lactam	Ampicillin-Sulbactam (f)	SXA	0.5/8-32/8	0.5/8-32/8	N/A
Beta-lactam	Ampicillin-Sulbactam	SAM	1/0.5-32/16	0.5/0.25-32/16	N/A
Aminoglycoside	Arbekacin	ARB	0.25-16	0.25-16	N/A
Macrolides Lincos-amides Streptogramins	Azithromycin	AZM	N/A	0.06-8	N/A

Antimicrobial			Available Concentrations		
Class	Drug Name ↓	Drug Code	Gram Negative	Gram Positive	Strep
			Range	Range	Range
Beta-lactam	Aztreonam	ATM	0.5-64	N/A	N/A
Beta-lactam	Cefaclor	CEC	N/A	0.5-32	N/A
Beta-lactam	Cefazolin	CZ	0.5-32	0.5-32	N/A
Beta-lactam	Cefdinir	CDR	0.12-4	0.12-4	N/A
Beta-lactam	Cefditoren	CDN	0.125-8	0.125-8	N/A
Beta-lactam	Cefepime	FEP	0.5-64	1-64	0.06-4
Beta-lactam	Cefetamet-pivoxil	CAT	0.25-16	N/A	N/A
Beta-lactam	Cefixime	CFM	0.125-8	N/A	N/A
Beta-lactam	Cefmetazole	CMZ	2-64	1-64	N/A
Beta-lactam	Cefoperazone	CFP	0.5-64	1-64	N/A
Beta-lactam	Cefoperazone-Sulbactam	SCP	0.5/8-64/8	N/A	N/A
Beta-lactam	Cefotaxime	CTX	0.5-64	1-64	0.06-4
Beta-lactam	Cefotetan	CTT	2-64	1-64	N/A
Beta-lactam	Cefotiam	CFT	0.5-64	0.5-64	N/A
Beta-lactam	Cefoxitin	FOX	0.5-64	1-64	N/A
Beta-lactam	Cefozopran	CFZ	N/A	1-64	N/A
Beta-lactam	Cefpirome	CPO	0.5-64	0.5-64	N/A
Beta-lactam	Cefpodoxime	CPD	0.12-8	0.5-8	N/A
Beta-lactam	Cefsulodin	CFS	1-64	N/A	N/A
Beta-lactam	Ceftazidime	CAZ	0.5-64	1-64	N/A
Beta-lactam	Ceftibuten	CTB	0.5-32	N/A	N/A
Beta-lactam	Ceftizoxime	ZOX	0.5-64	1-64	N/A
Beta-lactam	Ceftriaxone	CRO	0.5-64	1-64	0.06-4
Beta-lactam	Cefuroxime	CXM	1-64	1-64	0.12-4
Beta-lactam	Cephalexin	CN	1-64	0.5-64	N/A
Beta-lactam	Cephalothin	CF	1-64	0.5-64	N/A
Phenicol	Chloramphenicol	C	1-32	1-32	1-32

Antimicrobial			Available Concentrations		
Class	Drug Name ↓	Drug Code	Gram Negative	Gram Positive	Strep
			Range	Range	Range
Quinolone	Ciprofloxacin	CIP	0.12-4	0.12-4	N/A
Macrolides Lincos-amides Streptogramins	Clarithromycin	CLR	N/A	0.06-8	N/A
Macrolide Lincosamide	Clindamycin	CC	N/A	0.12-8	0.03-4
Cyclic peptide	Colistin	CL	0.5-4	N/A	N/A
Cyclic lipopeptide	Daptomycin	DAP	N/A	0.125-32	0.03-16
Tetracycline	Doxycycline	D	N/A	0.25-16	N/A
Beta-lactam	Ertapenem	ETP	0.25-32	0.25-32	0.06-4
Macrolide Lincosamide	Erythromycin	E	N/A	0.125-8	0.015-4
Fosfomycin	Fosfomycin	FF	8-256	8-256	N/A
Fusidane	Fusidic Acid	FA	N/A	0.5-32	N/A
Quinolone	Garenoxacin	GRN	0.125-16	0.25-8	0.03-4
Quinolone	Gatifloxacin	GAT	0.25-8	0.25-8	0.06-8
Quinolone	Gemifloxacin	GEM	0.125-8	0.125-2	0.06-8
Aminoglycoside	Gentamicin	GM	0.5-16	0.5-16	N/A
Aminoglycoside	Gentamicin-Syn	GMS	N/A	500	250-1000
Aminoglycoside	Gentamicin-Syn (SFM)	GMF	N/A	250	N/A
Beta-lactam	Imipenem	IPM	0.25-16	0.5-16	0.015-4
Aminoglycoside	Isepamycin	ISP	0.5-32	N/A	N/A
Aminoglycoside	Kanamycin	K	0.5-64	0.5-64	N/A
Aminoglycoside	Kanamycin-Syn	KS	N/A	250	250-1000
Quinolone	Levofloxacin	LVX	0.25-8	0.25-8	0.25-16
Macrolides Lincos-amides Streptogramin	Lincomycin	L	N/A	0.5-16	N/A
Oxazolidinone	Linezolid	LZD	N/A	0.25-32	0.25-16
Quinolone	Lomefloxacin	LOM	0.25-8	0.25-8	N/A
Beta-lactam	Mecillinam	MEC	0.5-32	N/A	N/A
Beta-lactam	Meropenem	MEM	0.25-16	0.5-16	0.03-2

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Class	Drug Name ↓	Drug Code	Range	Range	Range
Tetracycline	Minocycline	MI	0.5-16	0.5-32	N/A
Beta-lactam	Moxalactam	MOX	1-64	1-64	N/A
Quinolone	Moxifloxacin	MXF	0.12-8	0.12-8	0.06-8
Pseudomonic acid	Mupirocin	MUP	N/A	0.0625-8	N/A
Pseudomonic acid	Mupirocin High level	MUH	N/A	256	N/A
Quinolone	Nalidixic Acid	NA	1-32	N/A	N/A
Aminoglycoside	Netilmicin	NET	0.5-32	0.5-32	N/A
Nitrofurantoin	Nitrofurantoin	FM	8-512	16-512	N/A
Quinolone	Norfloxacin	NOR	0.25-16	0.25-16	N/A
Quinolone	Ofloxacin	OFX	0.25-8	0.25-8	0.5-16
Beta-lactam	Oxacillin	OX	N/A	0.06-4	N/A
Quinolone	Pefloxacin	PEF	0.25-8	0.25-8	N/A
Beta-lactam	Penicillin G	P	N/A	0.06-32	0.03-8
Beta-lactam	Piperacillin	PIP	0.5-128	1-128	N/A
Beta-lactam	Piperacillin-Tazobactam	TZP	0.5/4-128/4	1/4-128/4	N/A
Macrolides Lincosamides Streptogramin	Pristinamycin	PR	N/A	0.25-4	0.06-4
Macrolides Lincosamides Streptogramin	Quinupristin-Dalfopristin	SYN	N/A	0.5-4	0.12-8
Rifamycin	Rifampin	RA	N/A	0.25-32	N/A
Aminoglycoside	Streptomycin-Syn	STS	N/A	1000	250-1000
Aminoglycoside	Streptomycin-Syn (SFM)	STF	N/A	250	N/A
Glycopeptide	Teicoplanin	TEC	N/A	0.5-32	1-32
Ketolide	Telithromycin	TEL	N/A	0.03125-8	0.06-4
Beta-lactam	Temocillin	TEM	2 – 32	N/A	N/A
Tetracyclinee	Tetracycline	TE	0.5-16	0.5-16	0.06-16

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Class	Drug Name ↓	Drug Code	Range	Range	Range
Beta-lactam	Ticarcillin	TIC	1-128	1-128	N/A
Beta-lactam	Ticarcillin-Clavulanate	TIM	1/2-128/2	1/2-128/2	N/A
Glycylcycline	Tigecycline	TGC	0.25-16	N/A	N/A
Aminoglycoside	Tobramycin	NN	0.12-16	1 -16	N/A
Folate Antagonist	Trimeth-Sulfa (DIN)	STG	0.4/7.6-12.8/243.2	0.4/7.6-12.8/243.2	N/A
Folate Antagonist	Trimethoprim	TMP	0.5-16	0.5-16	N/A
Folate Antagonist	Trimethoprim-Sulfamethoxazole	SXT	0.5/9.5-16/304	0.5/9.5-16/304	0.06/1.19-16/304
Glycopeptide	Vancomycin	VA	N/A	0.5-32	0.06-32

BOLD = Different Range for Gram Negative and Gram Positive

In general, the Phoenix System provides a MIC for all organisms at any of the concentrations defined on a specific panel. For certain drug and organism combinations as indicated in the following tables, a specific minimum/maximum MIC is reported even if there is a lower/higher concentration on the panel.

Antimicrobial	Organism	Min MIC (µg/mL)	Max MIC (µg/mL)
Gram Negative			
Amikacin	<i>Morganella morganii</i>	2	
Amoxicillin-Clavulanic Acid (fixed clav. acid conc.)	<i>Proteus mirabilis</i>	1/2	
Amoxicillin-Clavulanic Acid (variable clav. acid conc.)	<i>Proteus mirabilis</i>	2/1	
Arbekacin	<i>Enterobacter aerogenes</i>	1	
	<i>Morganella morganii</i>	1	
	<i>Providencia</i> genus	1	
Cefdinir	<i>Proteus mirabilis</i>	0.25	
Cefixime	<i>Proteus mirabilis</i>	0.5	
	<i>Providencia</i> genus	0.5	
Cefotiam	<i>Citrobacter freundii</i>	2	
	<i>Proteus mirabilis</i>	4	
	<i>Providencia</i> genus	4	
Ceftriaxone	<i>Morganella morganii</i>	2	
Fosfomycin w/G6P	<i>Pseudomonas aeruginosa</i>	16	64
Fosfomycin w/G6P	<i>Serratia</i> genus	32	256
Garenoxacin	<i>Stenotrophomonas maltophilia</i>	1	
Imipenem	<i>Acinetobacter baumannii</i>	1	
	<i>Enterobacter aerogenes</i>	1	
	<i>Enterobacter cloacae</i>	1	
	<i>Morganella morganii</i>	2	
	<i>Proteus mirabilis</i>	2	
	<i>Proteus mirabilis</i>	4	
Kanamycin	<i>Proteus mirabilis</i>	4	
Meropenem	<i>Proteus mirabilis</i>	0.5	
Ticarclillin	<i>Proteus mirabilis</i>	4	

Antimicrobial	Organism	Min MIC (µg/mL)	Max MIC (µg/mL)
Gram Positive			
Amikacin	<i>Staphylococcus epidermidis</i>	8	
	<i>Staphylococcus</i> other*	8	
Amoxicillin	<i>Staphylococcus aureus</i>	2	
	<i>Staphylococcus epidermidis</i>	2	
Ampicillin	<i>Staphylococcus aureus</i>		1
	<i>Enterococcus</i> other**	1	
Ampicillin-Sulbactam (f) (SXA)	<i>Staphylococcus epidermidis</i>	1	
Ampicillin-Sulbactam (SAM)	<i>Enterococcus</i> other**	2	
Cefotetan	<i>Staphylococcus epidermidis</i>	8	
Ceftriaxone	<i>Staphylococcus aureus</i>	2	
Cephalexin	<i>Staphylococcus epidermidis</i>	1	
Daptomycin	<i>Enterococcus casseliflavus/</i> <i>gallinarum</i>	2	
	<i>Enterococcus</i> other**	1	8
Netilmicin	<i>Staphylococcus epidermidis</i>	2	
Norfloxacin	<i>Enterococcus</i> other**	2	
Penicillin	<i>Staphylococcus aureus</i>		1
	<i>Staphylococcus epidermidis</i>		1
	<i>Enterococcus</i> other**	2	
Piperacillin	<i>Staphylococcus epidermidis</i>	2	
Teicoplanin	<i>Staphylococcus epidermidis</i>	2	
Ticarcillin	<i>Staphylococcus epidermidis</i>	4	
	<i>Staphylococcus</i> other*	4	
Ticarcillin/Clavulanate	<i>Staphylococcus aureus</i>	4/2	
Trovafloxacin	<i>Enterococcus faecium</i>	1	
Streptococci			
Daptomycin	Viridans group	0.25	

* Coagulase negative *Staphylococcus* other than *S. epidermidis*

** *Enterococcus* sp. other than *E. faecalis* and *E. faecium*

F.2 List of Reagents and Principles Employed in the Phoenix System

F.2.1 Gram Negative

SUBSTRATE NAME	CODE	PRINCIPLE
L-PHENYLALANINE-AMC	A_LPHET	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
4MU-N-ACETYL-BD-GLUCOSAMINIDE	M_NAG	
L-GLUTAMIC ACID-AMC	A_LGTA	
L-TRYPTOPHAN-AMC	A_LTRY	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
L-PROLINE-AMC	A_LPROB	
L-ARGININE-AMC	A_LARGH	
ARGININE-ARGININE-AMC	A_ARARR	
GLYCINE-AMC	A_GLYB	
L-LEUCINE-AMC	A_LLEUH	
LYSINE-ALANINE-AMC	A_LYALD	
GLUTARYL-GLYCINE-ARGININE-AMC	A_GUGAH	
GLYCINE-PROLINE-AMC	A_GLPRB	
COLISTIN	C_CLST	Resistance to the antimicrobial agent results in a reduction of the resazurin based indicator.
POLYMYXIN B	C_PXB	
D-MANNITOL	C_DMNT	Utilization of a carbon source results in a reduction of the resazurin based indicator.
CITRATE	C_CIT	
ACETATE	C_ACT	
ADONITOL	C_ADO	
MALONATE	C_MLO	
ALPHA-KETOGLUTARIC ACID	C_KGA	
TIGLIC ACID	C_TIG	
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
L-PROLINE-NA	N_LPROT	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroaniline.
GAMMA-L-GLUTAMYL-NA	N_LGGH	
BIS (PNP) PHOSPHATE	P_BPHO	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol.
PNP-BD-GLUCOSIDE	P_BDGLU	

BETA-ALLOSE	R BALL	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
N-ACETYL-GALACTOSAMINE	R_NGA	
N-ACETYL-GLUCOSAMINE	R_NGU	
SORBITOL	R_DSBT	
SUCROSE	R_DSUC	
GALACTURONIC ACID	R_GRA	
MALTULOSE	R_MTU	
L-RHAMNOSE	R_LRHA	
BETA-GENTIOBIOSE	R_BGEN	
DEXTROSE	R_DEX	
D-GALACTOSE	R_DGAL	
D-FRUCTOSE	R_DFRU	
D-GLUCONIC ACID	R_DGUA	
D-MELIBIOSE	R_DMLB	
L-ARABINOSE	R_LARA	
METHYL-B-GLUCOSIDE	R_MBGU	
ORNITHINE	S_ORN	Utilization of ornithine results in pH rise and change in fluorescent indicator.
UREA	S_URE	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.

F.2.2 Gram Positive

SUBSTRATE NAME	CODE	PRINCIPLE
4MU-BD-CELLOBIOSIDE	M_BDCEL	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
L-ALANINE-AMC	A_LALT	
4MU-BD-GLUCOSIDE	M_BDGLU	
L-PROLINE-AMC	A_LPROB	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
L-PHENYLALANINE-AMC	A_LPHET	
L-TRYPTOPHAN-AMC	A_LTRY	
4MU-PHOSPHATE	M_PHOS	
METHIONINE-AMC	A_META	
4MU-AD-GLUCOSIDE	P_ADGLU	
ARGININE-ARGININE-AMC	A_ARARR	
GLYCINE-PROLINE-AMC	A_GLPRB	
4MU-BD-GLUCURONIDE	M_BDGLC	
L-LEUCINE-AMC	A_LLEUH	
4MU-N-ACETYL-BD-GLUCOSAMINIDE	M_NAG	
L-ARGININE-AMC	A_LARGH	
4MU-PHOSPHATE (with Trehalose)	M_PHOT	
L-HISTIDINE-AMC	A_LHIST	
L-ISOLEUCINE-AMC	A_LISO	
4MU-BD-GALACTOSIDE	M_BDGAL	

COLISTIN	C_CLST	}	Resistance to the antimicrobial agent results in a reduction of the resazurin based indicator.
POLYMYXIN B	C_PXB		
D-GLUCONIC ACID	C_DGUA	}	Utilization of a carbon source results in a reduction of the resazurin based indicator.
3-METHYL GLUTARIC ACID	C_3MGA		
D-FRUCTOSE	C_DFRU		
IMINODIACETIC ACID	C_IMN		
ALPHA-KETOGLUTARIC ACID	C_KGA		
D-MANNITOL	C_DMNT		
3-METHYLADIPIC ACID	C_MAA		
THYMIDINE	C_THY		
FLUORESCENT POSITIVE CONTROL	FLR_CTL	}	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL		
ALANINE-ALANINE-PNA	N_ALALH	}	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroaniline.
L-PROLINE-PNA	N_LPROT		
VALINE-ALANINE-PNA	N_VAALA		
PNP-AD-GLUCOSIDE	P_PAGLU	}	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol.
PNP-PHOSPHATE	P_PHOL		
BETA-GENTIOBIOSE	R_BGEN	}	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
D-SUCROSE	R_DSUC		
MALTOTRIOSE	R_MTT		
N-ACETYL-GLUCOSAMINE	R_NGU		
D-TREHALOSE	R_DTRE		
D-TAGATOSE	R_DTAG		
MALTOSE	R_MAL		
DEXTROSE	R_DEX		
METHYL-α-D-GLUCOSIDE	R_MGP		
UREA	S_URE	}	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.
ESCULIN	T_ESC	}	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.
NITROCEFIN	L_NCF	}	Enzymatic hydrolysis of the b-Lactam ring results in a color change

F.2.3 Streptococci Panel

SUBSTRATE NAME	CODE	PRINCIPLE
AMYGDALIN	R_AMY	Utilization of carbohydrate results in lower pH and change in indicator (Phenol red).
D-GALACTOSE	R_DGAL	
D-MANNITOL	R_DMTL	
D-RAFFINOSE	R_DRAF	
D-SORBITOL	R_DSBT	
D-TREHALOSE	R_DTRE	
DEXTRIN	R_DNX	
N-ACETYL-GLUCOSAMINE	R_NGU	
PHENYL GLUCOSIDE	R_PHG	
SALICIN	R_SAL	
ONP-BD-GLUCOSIDE	O_BOGLU	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol.
PNP-AD-GALACTOSIDE	P_ADGAL	
PNP-BD-CELLOBIOSIDE	P_CELB	
PNP-BD-GALACTOSIDE	P_GALB	
PNP-AD-GLUCOSIDE	P_PAGLU	
PNP-PHOSPHATE	P_PHOL	
ALANINE-ALANINE-PNA	N_ALALH	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroanilide.
VALINE-ALANINE-PNA	N_VAALA	
L-LYSINE-PNA	N_LLYSB	
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
THYMIDINE	THY	Utilization of a carbon source resulting in a reduction of the indicator (Resazurin based).
PULLULAN	PUL	
D-TREHALOSE	TRL	
D-LACTOSE	DLAC	

LYSINE-AMC	A_LYSA	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
SERINE-TYROSINE-AMC	A_SETY	
L-CITRULLINE-AMC	A_LCTU	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
ISOLEUCINE-AMC	A_LISO	
L-TRYPTOPHAN-AMC	A_LTRY	
L-VALINE-AMC	A_LVAL	
ARGININE-ARGININE-AMC	A_ARARR	
LYSINE-ALANINE-AMC	A_LYALD	
ASPARAGINE-AMC	A_APGT	
L-ARGININE-AMC	A_LARGH	
L-HISTIDINE-AMC	A_LHIST	
ALANINE-AFC	Z_ALFT	
4MU-BD-CELLOBIOSIDE	M_BDCEL	
4MU-BD-GLUCOSIDE	M_BDGLU	
4MU-PHOSPHATE	M_PHOS	
4MU-AD-GLUCOSIDE	M_ADGLU	
4MU-BD-GLUCURONIDE	M_BDGLC	
4MU-N-ACETYL-BD-GLUCOSAMINE	M_NAG	
4MU-PHOSPHATE (with trehalose)	M_PHOT	
4MU-BD-GALACTOSIDE	M_BDGAL	
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.

F.2.4 Yeast Panel

SUBSTRATE NAME	CODE	PRINCIPLE
PNP-BD-GLUCOSIDE	P_BDGLU	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow p-nitrophenol.
PNP-AD-GLUCOSIDE	P_PAGLU	
ONP-BD-GLUCOSIDE	O_BOGLU	
L-SORBOSE	C_LSBO	Utilization of a carbon source results in a reduction of the resazurin based indicator.
DEXTROSE	C_DEX	
D-MANNITOL	C_DMNT	
D-SUCROSE	C_DSUC	
METHYL-AD-GLUCOPYRANOSIDE	C_MGP	
N-ACETYL-BD-GLUCOSAMINIDE	C_NAG	
DEXTROSE	R_DEX	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
D-FRUCTOSE	R_DFRU	
D-GALACTOSE	R_DGAL	
SUCROSE	R_DSUC	
D-TREHALOSE	R_DTRE	
MALTOTRIOSE	R_MTT	
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.
FLUORESCENT NEGATIVE CONTROL	Z_FTST	Control to check for fluorescent interference.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
GAMMA-L-GLUTAMYL-NA	N_LGGH	Enzymatic hydrolysis of the colorless amide substrate releases yellow p-nitroaniline.
L-PROLINE-PNA	N_LPROT	

ASPARAGINE-AMC	A_APGT	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
L-ARGININE-AMC	A_LARGH	
L-GLUTAMINE-AMC	A_LGLNB	
L-TYROSINE-AMC	A_LTYO	
L-HISTIDINE-AMC	A_LHIST	
ORNITHINE-AMC	A_ORN	
THREONINE-AMC	A_THR	
HYDROXY PROLINE-AMC	A_LHYP	
4MU-N-ACETYL-BD-GLUCOSAMINE	M_NAG	
4MU-AD-GLUCOSIDE	M_ADGLU	
4MU-PHOSPHATE	M_PHOS	
LYSINE-ALANINE-AMC	A_LYALD	
GLYCINE-ARGININE-AMC	A_GLARH	
ALANINE-AFC	Z_ALFT	
GLYCINE-AMC	A_GLYB	
L-CITRULLINE-AMC	A_LCTU	
L-GLUTAMIC ACID-AMC	A_LGTA	
L-VALINE-AMC	A_LVAL	
L-ALANINE-AMC	A_LALT	
L-PROLINE-AMC	A_LPROB	
L-TRYPTOPHAN-AMC	A_LTRY	
H-B-ALANINE-AMC	A_HBALT	
4MU-BD-CELLOBIOSIDE	M_BDCEL	
4MU-BD-GLUCOSIDE	M_BDGLU	
GLYCINE-PROLINE-AMC	A_GLPRB	
LYSINE-PROLINE-AMC	A_LYPRA	
BENZYL-L-CYSTEINE-AMC	A_BZLCY	
AMINO ACID	S_GTN	Utilization of the amino acid results in a change in fluorescence.
UREA	S_URE	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.

F.3 Taxa for ID/AST Determination

F.3.1 Gram Negative (0.5 McFarland)

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST ²	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST ²
<i>Achromobacter denitrificans</i>	AST	<i>Aeromonas schubertii</i>	ID/AST
<i>Achromobacter piechaudii</i>	AST	<i>Aeromonas sobria</i>	ID/AST
<i>Achromobacter</i> species	ID/AST	<i>Aeromonas trota</i>	AST
<i>Achromobacter xylosoxidans</i> ssp <i>xylosoxidans</i>	AST	<i>Aeromonas veronii</i>	ID/AST
<i>Acinetobacter baumannii</i> ²	ID/AST	<i>Alcaligenes faecalis</i>	ID/AST
<i>Acinetobacter baumannii</i> / <i>calcoaceticus</i> complex ²	ID/AST	<i>Bergeyella zoohelcum</i>	ID
<i>Acinetobacter calcoaceticus</i>	AST	<i>Bordetella bronchiseptica</i>	ID
<i>Acinetobacter haemolyticus</i>	ID/AST	<i>Brevundimonas diminuta</i>	ID/AST
<i>Acinetobacter johnsonii</i>	AST	<i>Brevundimonas vesicularis</i>	ID/AST
<i>Acinetobacter junii</i>	AST	<i>Burkholderia caryophylli</i>	AST
<i>Acinetobacter lwoffii</i>	ID/AST	<i>Burkholderia cepacia</i>	ID/AST
<i>Acinetobacter radioresistens</i>	AST	<i>Burkholderia gladioli</i>	ID/AST
<i>Acinetobacter</i> species	ID/AST	<i>Burkholderia glathei</i>	AST
<i>Actinobacillus lignieresii</i>	ID	<i>Burkholderia graminis</i>	AST
<i>Actinobacillus suis</i>	ID	<i>Burkholderia multivorans</i>	AST
<i>Actinobacillus ureae</i>	ID	<i>Burkholderia phenazinium</i>	AST
<i>Aeromonas allosaccharophila</i>	AST	<i>Burkholderia pyrrocinia</i>	AST
<i>Aeromonas caviae</i>	ID/AST	<i>Cardiobacterium hominis</i>	ID
<i>Aeromonas eucrenophila</i>	AST	CDC group EF-4a	ID
<i>Aeromonas hydrophila</i>	ID/AST	CDC group EF-4b	ID
<i>Aeromonas jandaei</i>	AST	CDC group EO-2	ID
<i>Aeromonas media</i>	AST	CDC group Vb-3	ID
<i>Aeromonas salmonicida</i>	AST	<i>Cedecea davisae</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>achromogenes</i>	AST	<i>Cedecea lapagei</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>masoucida</i>	ID/AST	<i>Cedecea neteri</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>pectinolytica</i>	AST	<i>Cedecea</i> species 3	AST
<i>Aeromonas salmonicida</i> ssp <i>salmonicida</i>	ID/AST	<i>Cedecea</i> species 5	AST
<i>Aeromonas salmonicida</i> ssp <i>smithia</i>	ID/AST	<i>Chromobacterium violaceum</i>	ID
		<i>Chryseobacterium gleum</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).
¹ Not all species encountered during clinical performance evaluations
² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Chryseobacterium indologenes</i>	ID/AST	<i>Enterobacter nimipressuralis</i>	AST
<i>Chryseobacterium scopthalmum</i>	AST	<i>Enterobacter sakazakii</i>	ID/AST
<i>Citrobacter amalonaticus</i>	ID/AST	<i>Escherichia blattae</i>	AST
<i>Citrobacter braakii</i>	ID/AST	<i>Escherichia coli</i> ²	ID/AST
<i>Citrobacter farmeri</i>	ID/AST	<i>Escherichia fergusonii</i>	ID/AST
<i>Citrobacter freundii</i>	ID/AST	<i>Escherichia hermannii</i>	ID/AST
<i>Citrobacter gillenii</i>	AST	<i>Escherichia vulneris</i>	ID/AST
<i>Citrobacter koseri</i>	ID/AST	<i>Ewingella americana</i>	ID
<i>Citrobacter murlinae</i>	AST	<i>Grimontia hollisae</i>	ID
<i>Citrobacter rodentium</i>	AST	<i>Hafnia alvei</i>	ID/AST
<i>Citrobacter sedlakii</i>	ID/AST	<i>Kingella denitrificans</i>	ID
<i>Citrobacter werkmanii</i>	ID/AST	<i>Kingella kingae</i>	ID
<i>Citrobacter youngae</i>	ID/AST	<i>Klebsiella granulomatis</i>	AST
<i>Comamonas terrigena</i>	ID	<i>Klebsiella oxytoca</i> ²	ID/AST
<i>Comamonas testosteroni</i>	ID	<i>Klebsiella pneumoniae</i> ssp <i>ozaenae</i>	ID/AST
<i>Cupriavidus gilardii</i>	AST	<i>Klebsiella pneumoniae</i> ssp <i>pneumoniae</i> ²	ID/AST
<i>Cupriavidus pauculus</i>	ID/AST	<i>Klebsiella pneumoniae</i> ssp <i>rhinoscleromatis</i>	ID/AST
<i>Delftia acidovorans</i>	ID/AST	<i>Kluyvera ascorbata</i>	ID/AST
<i>Edwardsiella hoshinae</i>	ID/AST	<i>Kluyvera cryocrescens</i>	ID/AST
<i>Edwardsiella ictaluri</i>	ID/AST	<i>Kluyvera georgiana</i>	AST
<i>Edwardsiella tarda</i>	ID/AST	<i>Kluyvera intermedia</i>	ID/AST
<i>Eikenella corrodens</i>	ID	<i>Leclercia adecarboxylata</i>	ID/AST
<i>Elizabethkingia meningoseptica</i>	ID/AST	<i>Leminorella grimontii</i>	ID
<i>Empedobacter brevis</i>	ID	<i>Leminorella richardii</i>	ID
<i>Enterobacter aerogenes</i> ²	ID/AST	<i>Mannheimia haemolytica</i>	ID
<i>Enterobacter amnigenus</i>	AST	<i>Methylobacterium extorquens</i>	ID
<i>Enterobacter amnigenus</i> biogroup 1	ID/AST	<i>Moellerella wisconsensis</i>	ID/AST
<i>Enterobacter amnigenus</i> biogroup 2	ID/AST	<i>Moraxella (Branhamella)</i> <i>catarrhalis</i>	ID
<i>Enterobacter asburiae</i>	ID/AST	<i>Moraxella species</i>	ID
<i>Enterobacter cancerogenus</i>	ID/AST	<i>Morganella morganii</i> ²	ID/AST
<i>Enterobacter cloacae</i> ²	ID/AST	<i>Myroides odoratus/odoratimimus</i>	ID/AST
<i>Enterobacter cloacae</i> ssp <i>dissolvens</i>	AST	<i>Ochrobactrum anthropi</i>	ID/AST
<i>Enterobacter cowanii</i>	AST	<i>Oligella ureolytica</i>	ID
<i>Enterobacter gergoviae</i>	ID/AST	<i>Oligella urethralis</i>	ID
<i>Enterobacter hormaechei</i>	ID/AST	<i>Pantoea agglomerans</i>	ID/AST
<i>Enterobacter kobei</i>	AST	<i>Pantoea ananatis</i>	AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Pantoea dispersa</i>	AST	<i>Ralstonia species</i>	AST
<i>Pantoea stewartii</i> ssp <i>indologenes</i>	AST	<i>Raoultella ornithinolytica</i>	ID/AST
<i>Pantoea stewartii</i> ssp <i>stewartii</i>	AST	<i>Raoultella planticola</i>	AST
<i>Pasteurella aerogenes</i>	ID	<i>Raoultella terrigena</i>	AST
<i>Pasteurella multocida</i>	ID	<i>Rhizobium radiobacter</i>	ID
<i>Pasteurella pneumotropica</i>	ID	<i>Salmonella aberdeen</i>	AST
<i>Photobacterium damsela</i>	ID	<i>Salmonella abortus-equi</i>	AST
<i>Plesiomonas shigelloides</i>	ID	<i>Salmonella adelaide</i>	AST
<i>Pragia fontium</i>	ID	<i>Salmonella aderike</i>	AST
<i>Proteus hauseri</i>	AST	<i>Salmonella agona</i>	AST
<i>Proteus mirabilis</i> ²	ID/AST	<i>Salmonella alachua</i>	AST
<i>Proteus myxofaciens</i>	AST	<i>Salmonella anatum</i>	AST
<i>Proteus penneri</i>	ID/AST	<i>Salmonella arizonae</i>	AST
<i>Proteus vulgaris</i>	ID/AST	<i>Salmonella avana</i>	AST
<i>Providencia alcalifaciens</i>	ID/AST	<i>Salmonella bahrenfeld</i>	AST
<i>Providencia heimbachae</i>	AST	<i>Salmonella blockley</i>	AST
<i>Providencia rettgeri</i>	ID/AST	<i>Salmonella bongori</i>	AST
<i>Providencia rustigianii</i>	ID/AST	<i>Salmonella braenderup</i>	AST
<i>Providencia stuartii</i>	ID/AST	<i>Salmonella bredeney</i>	AST
<i>Pseudomonas aeruginosa</i> ²	ID/AST	<i>Salmonella bunn</i>	AST
<i>Pseudomonas alcaligenes</i>	AST	<i>Salmonella californica</i>	AST
<i>Pseudomonas fluorescens</i>	ID/AST	<i>Salmonella carrau</i>	AST
<i>Pseudomonas luteola</i>	ID/AST	<i>Salmonella cerro</i>	AST
<i>Pseudomonas mendocina</i>	ID/AST	<i>Salmonella champaign</i>	AST
<i>Pseudomonas monteilii</i>	AST	<i>Salmonella chittagong</i>	AST
<i>Pseudomonas oryziatrans</i>	ID/AST	<i>Salmonella choleraesuis</i>	AST
<i>Pseudomonas pertucinogena</i>	AST	<i>Salmonella choleraesuis</i> ssp <i>arizonae</i>	ID/AST
<i>Pseudomonas pseudoalcaligenes</i>	ID/AST	<i>Salmonella choleraesuis</i> ssp <i>choleraesuis</i>	ID/AST
<i>Pseudomonas putida</i>	ID/AST	<i>Salmonella choleraesuis</i> ssp <i>diarizonae</i>	AST
<i>Pseudomonas species</i>	ID/AST	<i>Salmonella choleraesuis</i> ssp <i>houtenae</i>	AST
<i>Pseudomonas stutzeri</i>	ID/AST	<i>Salmonella choleraesuis</i> ssp <i>indica</i>	AST
<i>Pseudomonas veronii</i>	AST	<i>Salmonella choleraesuis</i> ssp <i>salamae</i>	AST
<i>Rahnella aquatilis</i>	ID	<i>Salmonella cubana</i>	AST
<i>Ralstonia pickettii</i>	ID/AST	<i>Salmonella dakar</i>	AST
<i>Ralstonia solanacearum</i>	AST	<i>Salmonella daressalaam</i>	AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Salmonella derby</i>	AST	<i>Salmonella montevideo</i>	AST
<i>Salmonella dessau</i>	AST	<i>Salmonella muenchen</i>	AST
<i>Salmonella DT</i>	AST	<i>Salmonella muenster</i>	AST
<i>Salmonella dublin</i>	AST	<i>Salmonella newington</i>	AST
<i>Salmonella duesseldorf</i>	AST	<i>Salmonella newport</i>	AST
<i>Salmonella enteritidis</i>	AST	<i>Salmonella nottingham</i>	AST
<i>Salmonella fresno</i>	AST	<i>Salmonella ohio</i>	AST
<i>Salmonella gallinarum</i>	ID/AST	<i>Salmonella onderstepoort</i>	AST
<i>Salmonella gaminara</i>	AST	<i>Salmonella oranienburg</i>	AST
<i>Salmonella give</i>	AST	<i>Salmonella panama</i>	AST
<i>Salmonella haardt</i>	AST	<i>Salmonella paratyphi A</i>	ID/AST
<i>Salmonella hadar</i>	AST	<i>Salmonella paratyphi B</i>	AST
<i>Salmonella hamburg</i>	AST	<i>Salmonella poona</i>	AST
<i>Salmonella hartford</i>	AST	<i>Salmonella pullorum</i>	ID/AST
<i>Salmonella heidelberg</i>	AST	<i>Salmonella quinhon</i>	AST
<i>Salmonella illinois</i>	AST	<i>Salmonella rubislaw</i>	AST
<i>Salmonella infantis</i>	AST	<i>Salmonella saintpaul</i>	AST
<i>Salmonella inverness</i>	AST	<i>Salmonella schwarzengrund</i>	AST
<i>Salmonella java</i>	AST	<i>Salmonella senftenberg</i>	AST
<i>Salmonella javiana</i>	AST	<i>Salmonella species</i>	ID/AST
<i>Salmonella kentucky</i>	AST	<i>Salmonella tallahassee</i>	AST
<i>Salmonella kirkee</i>	AST	<i>Salmonella thompson</i>	AST
<i>Salmonella kunduchi</i>	AST	<i>Salmonella typhi</i>	ID/AST
<i>Salmonella kvittingfoss</i>	AST	<i>Salmonella typhimurium</i>	AST
<i>Salmonella lansing</i>	AST	<i>Salmonella virginia</i>	AST
<i>Salmonella litchfield</i>	AST	<i>Salmonella westerstede</i>	AST
<i>Salmonella liverpool</i>	AST	<i>Salmonella worthington</i>	AST
<i>Salmonella london</i>	AST	<i>Serratia entomophila</i>	AST
<i>Salmonella luciana</i>	AST	<i>Serratia ficaria</i>	ID/AST
<i>Salmonella manhattan</i>	AST	<i>Serratia fonticola</i>	ID/AST
<i>Salmonella mbandaka</i>	AST	<i>Serratia grimesii</i>	AST
<i>Salmonella meleagridis</i>	AST	<i>Serratia liquefaciens</i>	ID/AST
<i>Salmonella memphis</i>	AST	<i>Serratia marcescens</i> ²	ID/AST
<i>Salmonella michigan</i>	AST	<i>Serratia odorifera</i>	AST
<i>Salmonella minneapolis</i>	AST	<i>Serratia odorifera 1</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Serratia odorifera</i> 2	ID/AST	<i>Vibrio alginolyticus</i>	ID
<i>Serratia plymuthica</i>	ID/AST	<i>Vibrio cholerae</i>	ID
<i>Serratia proteamaculans</i> ssp <i>proteamaculans</i>	AST	<i>Vibrio fluvialis</i>	ID
<i>Serratia proteamaculans</i> ssp <i>quinovora</i>	AST	<i>Vibrio metschnikovii</i>	ID
<i>Serratia rubidaea</i>	ID/AST	<i>Vibrio mimicus</i>	ID
<i>Shewanella algae</i>	AST	<i>Vibrio parahaemolyticus</i>	ID
<i>Shewanella putrefaciens</i>	ID/AST	<i>Vibrio vulnificus</i>	ID
<i>Shigella boydii</i>	ID/AST	<i>Weeksella virosa</i>	ID
<i>Shigella dysenteriae</i>	ID/AST	<i>Yersinia aldovae</i>	AST
<i>Shigella flexneri</i>	ID/AST	<i>Yersinia bercovieri</i>	AST
<i>Shigella sonnei</i>	ID/AST	<i>Yersinia enterocolitica</i>	ID/AST
<i>Shigella</i> species	ID/AST	<i>Yersinia frederiksenii</i>	ID/AST
<i>Sphingobacterium multivorum</i>	ID/AST	<i>Yersinia intermedia</i>	ID/AST
<i>Sphingobacterium spiritivorum</i>	ID/AST	<i>Yersinia kristensenii</i>	ID/AST
<i>Sphingobacterium thalpophilum</i>	ID/AST	<i>Yersinia mollaretii</i>	AST
<i>Sphingomonas paucimobilis</i>	ID/AST	<i>Yersinia pseudotuberculosis</i>	ID/AST
<i>Stenotrophomonas africana</i>	AST	<i>Yersinia rohdei</i>	AST
<i>Stenotrophomonas maltophilia</i> ²	ID/AST	<i>Yersinia ruckeri</i>	ID/AST
<i>Suttonella indologenes</i>	ID	<i>Yokenella regensburgei</i>	ID
<i>Tatumella ptyseos</i>	ID		

F.3.2 Gram Negative (0.25 McFarland)

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Achromobacter</i> species	ID/AST	<i>Aeromonas caviae</i>	ID/AST
<i>Acinetobacter baumannii</i>		<i>Aeromonas hydrophila</i>	ID/AST
<i>calcoaceticus</i> complex	ID/AST	<i>Aeromonas salmonicida</i> ssp <i>masoucida</i>	ID/AST
<i>Acinetobacter haemolyticus</i>	ID/AST	<i>Aeromonas salmonicida</i> ssp <i>salmonicida</i>	ID/AST
<i>Acinetobacter lwoffii</i>	ID/AST	<i>Aeromonas salmonicida</i> ssp <i>smithia</i>	ID/AST
<i>Actinobacillus lignieresii</i>	ID	<i>Aeromonas schubertii</i>	ID/AST
<i>Actinobacillus suis</i>	ID	<i>Aeromonas sobria</i>	ID/AST
<i>Actinobacillus ureae</i>	ID	<i>Aeromonas veronii</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Alcaligenes faecalis</i>	ID/AST	<i>Eikenella corrodens</i>	ID
<i>Bergeyella zoohelcum</i>	ID	<i>Elizabethkingia meningoseptica</i>	ID/AST
<i>Bordetella bronchiseptica</i>	ID	<i>Empedobacter brevis</i>	ID
<i>Brevundimonas diminuta</i>	ID/AST	<i>Enterobacter aerogenes</i>	ID/AST
<i>Brevundimonas vesicularis</i>	ID/AST	<i>Enterobacter amnigenus</i> biogroup 1	ID/AST
<i>Burkholderia cepacia</i>	ID/AST	<i>Enterobacter amnigenus</i> biogroup 2	ID/AST
<i>Burkholderia gladioli</i>	ID/AST	<i>Enterobacter asburiae</i>	ID/AST
<i>Cardiobacterium hominis</i>	ID	<i>Enterobacter cancerogenus</i>	ID/AST
CDC group EF-4a	ID	<i>Enterobacter cloacae</i>	ID/AST
CDC group EF-4b	ID	<i>Enterobacter gergoviae</i>	ID/AST
CDC group EO-2	ID	<i>Enterobacter hormaechei</i>	ID/AST
CDC group Vb-3	ID	<i>Enterobacter sakazakii</i>	ID/AST
<i>Cedecea davisae</i>	ID/AST	<i>Escherichia coli</i>	ID/AST
<i>Cedecea lapagei</i>	ID/AST	<i>Escherichia fergusonii</i>	ID/AST
<i>Cedecea neteri</i>	ID/AST	<i>Escherichia hermannii</i>	ID/AST
<i>Chromobacterium violaceum</i>	ID	<i>Escherichia vulneris</i>	ID/AST
<i>Chryseobacterium gleum</i>	ID/AST	<i>Ewingella americana</i>	ID
<i>Chryseobacterium indologenes</i>	ID/AST	<i>Grimontia hollisae</i>	ID
<i>Citrobacter amalonaticus</i>	ID/AST	<i>Hafnia alvei</i>	ID/AST
<i>Citrobacter braakii</i>	ID/AST	<i>Klebsiella oxytoca</i>	ID/AST
<i>Citrobacter farmeri</i>	ID/AST	<i>Klebsiella pneumoniae</i> ssp ozaenae	ID/AST
<i>Citrobacter freundii</i>	ID/AST	<i>Klebsiella pneumoniae</i> ssp pneumoniae	ID/AST
<i>Citrobacter koseri</i>	ID/AST	<i>Klebsiella pneumoniae</i> ssp rhinoscleromatis	ID/AST
<i>Citrobacter sedlakii</i>	ID/AST	<i>Kluyvera ascorbata</i>	ID/AST
<i>Citrobacter werkmanii</i>	ID/AST	<i>Kluyvera cryocrescens</i>	ID/AST
<i>Citrobacter youngae</i>	ID/AST	<i>Kluyvera intermedia</i>	ID/AST
<i>Comamonas terrigena</i>	ID	<i>Leclercia adecarboxylata</i>	ID/AST
<i>Comamonas testosteroni</i>	ID	<i>Leminorella grimontii</i>	ID
<i>Cupriavidus pauculus</i>	ID/AST	<i>Leminorella richardii</i>	ID
<i>Delftia acidovorans</i>	ID/AST	<i>Mannheimia haemolytica</i>	ID
<i>Edwardsiella hoshinae</i>	ID/AST	<i>Moellerella wisconsensis</i>	ID/AST
<i>Edwardsiella ictaluri</i>	ID/AST	<i>Morganella morganii</i>	ID/AST
<i>Edwardsiella tarda</i>	ID/AST	<i>Myroides odoratus/odoratimimus</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).
¹ Not all species encountered during clinical performance evaluations

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Ochrobactrum anthropi</i>	ID/AST	<i>Serratia ficaria</i>	ID/AST
<i>Oligella ureolytica</i>	ID	<i>Serratia fonticola</i>	ID/AST
<i>Oligella urethralis</i>	ID	<i>Serratia liquefaciens</i>	ID/AST
<i>Pantoea agglomerans</i>	ID/AST	<i>Serratia marcescens</i>	ID/AST
<i>Pasteurella aerogenes</i>	ID	<i>Serratia odorifera</i> 1	ID/AST
<i>Pasteurella multocida</i>	ID	<i>Serratia odorifera</i> 2	ID/AST
<i>Pasteurella pneumotropica</i>	ID	<i>Serratia plymuthica</i>	ID/AST
<i>Photobacterium damsela</i>	ID	<i>Serratia rubidaea</i>	ID/AST
<i>Plesiomonas shigelloides</i>	ID	<i>Shewanella putrefaciens</i>	ID/AST
<i>Pragia fontium</i>	ID	<i>Shigella boydii</i>	ID/AST
<i>Proteus mirabilis</i>	ID/AST	<i>Shigella dysenteriae</i>	ID/AST
<i>Proteus penneri</i>	ID/AST	<i>Shigella flexneri</i>	ID/AST
<i>Proteus vulgaris</i>	ID/AST	<i>Shigella sonnei</i>	ID/AST
<i>Providencia alcalifaciens</i>	ID/AST	<i>Sphingobacterium multivorum</i>	ID/AST
<i>Providencia rettgeri</i>	ID/AST	<i>Sphingobacterium spiritivorum</i>	ID/AST
<i>Providencia rustigianii</i>	ID/AST	<i>Sphingobacterium thalpophilum</i>	ID/AST
<i>Providencia stuartii</i>	ID/AST	<i>Sphingomonas paucimobilis</i>	ID/AST
<i>Pseudomonas aeruginosa</i>	ID/AST	<i>Stenotrophomonas maltophilia</i>	ID/AST
<i>Pseudomonas fluorescens</i>	ID/AST	<i>Suttonella indologenes</i>	ID
<i>Pseudomonas luteola</i>	ID/AST	<i>Tatumella ptyseos</i>	ID
<i>Pseudomonas mendocina</i>	ID/AST	<i>Vibrio alginolyticus</i>	ID
<i>Pseudomonas oryzae</i>	ID/AST	<i>Vibrio cholerae</i>	ID
<i>Pseudomonas putida</i>	ID/AST	<i>Vibrio fluvialis</i>	ID
<i>Pseudomonas stutzeri</i>	ID/AST	<i>Vibrio metschnikovii</i>	ID
<i>Rahnella aquatilis</i>	ID	<i>Vibrio mimicus</i>	ID
<i>Ralstonia pickettii</i>	ID/AST	<i>Vibrio parahaemolyticus</i>	ID
<i>Raoultella ornithinolytica</i>	ID/AST	<i>Vibrio vulnificus</i>	ID
<i>Rhizobium radiobacter</i>	ID	<i>Weeksella virosa</i>	ID
<i>Salmonella choleraesuis</i> ssp <i>arizonae</i>	ID/AST	<i>Yersinia enterocolitica</i>	ID/AST
<i>Salmonella choleraesuis</i> ssp <i>choleraesuis</i>	ID/AST	<i>Yersinia frederiksenii</i>	ID/AST
<i>Salmonella gallinarum</i>	ID/AST	<i>Yersinia intermedia</i>	ID/AST
<i>Salmonella paratyphi</i> A	ID/AST	<i>Yersinia kristensenii</i>	ID/AST
<i>Salmonella pullorum</i>	ID/AST	<i>Yersinia pseudotuberculosis</i>	ID/AST
<i>Salmonella species</i>	ID/AST	<i>Yersinia ruckeri</i>	ID/AST
<i>Salmonella typhi</i>	ID/AST	<i>Yokenella regensburgei</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

F.3.3 Gram Positive (0.5 McFarland)

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Aerococcus urinae</i>	ID	<i>Enterococcus casseliflavus</i>	ID/AST
<i>Aerococcus viridans</i>	ID	<i>Enterococcus cecorum</i>	AST
<i>Alloiococcus otitidis</i>	ID	<i>Enterococcus columbae</i>	AST
<i>Arcanobacterium haemolyticum</i>	ID	<i>Enterococcus dispar</i>	AST
<i>Arcanobacterium pyogenes</i>	ID	<i>Enterococcus durans</i>	ID/AST
<i>Bacillus cereus</i>	ID	<i>Enterococcus faecalis</i>	ID/AST
<i>Bacillus circulans</i>	ID	<i>Enterococcus faecium</i>	ID/AST
<i>Bacillus coagulans</i>	ID	<i>Enterococcus flavescens</i>	AST
<i>Bacillus licheniformis</i>	ID	<i>Enterococcus gallinarum</i>	ID/AST
<i>Bacillus megaterium</i>	ID	<i>Enterococcus gilvus</i>	AST
<i>Bacillus pumilus</i>	ID	<i>Enterococcus haemoperoxidus</i>	AST
<i>Bacillus subtilis</i>	ID	<i>Enterococcus hirae</i>	ID/AST
<i>Bacillus thuringiensis</i>	ID	<i>Enterococcus malodoratus</i>	AST
<i>Brevibacillus brevis</i>	ID	<i>Enterococcus moraviensis</i>	AST
<i>Brevibacterium species</i>	ID	<i>Enterococcus mundtii</i>	AST
<i>Cellulomonas turbata</i>	ID	<i>Enterococcus pallens</i>	AST
<i>Corynebacterium amycolatum</i>	ID	<i>Enterococcus pseudoavium</i>	AST
<i>Corynebacterium bovis</i>	ID	<i>Enterococcus raffinosus</i>	ID/AST
<i>Corynebacterium diphtheriae</i>	ID	<i>Enterococcus ratti</i>	AST
<i>Corynebacterium jeikeium</i>	ID	<i>Enterococcus saccharolyticus</i>	AST
<i>Corynebacterium kutscheri</i>	ID	<i>Enterococcus sulfureus</i>	AST
<i>Corynebacterium matruchotii</i>	ID	<i>Erysipelothrix rhusiopathiae</i>	ID
<i>Corynebacterium minutissimum</i>	ID	<i>Gardnerella vaginalis</i>	ID
<i>Corynebacterium propinquum</i>	ID	<i>Gemella haemolysans</i>	ID
<i>Corynebacterium pseudodiphtheriticum</i>	ID	<i>Gemella morbillorum</i>	ID
<i>Corynebacterium pseudotuberculosis</i>	ID	<i>Globicatella sanguinis</i>	ID
<i>Corynebacterium renale</i>	ID	<i>Helcococcus kunzii</i>	ID
<i>Corynebacterium striatum</i>	ID	<i>Kocuria kristinae</i>	ID
<i>Corynebacterium ulcerans</i>	ID	<i>Kocuria rosea</i>	ID
<i>Corynebacterium urealyticum</i>	ID	<i>Kocuria varians</i>	ID
<i>Corynebacterium xerosis</i>	ID	<i>Kytococcus sedentarius</i>	ID
<i>Dermabacter hominis</i>	ID	<i>Lactococcus garvieae</i>	ID
<i>Dermacoccus nishinomiyaensis</i>	ID	<i>Lactococcus lactis ssp cremoris</i>	ID
<i>Enterococcus asini</i>	AST	<i>Lactococcus lactis ssp hordniae</i>	ID
<i>Enterococcus avium</i>	ID/AST	<i>Lactococcus lactis ssp lactis</i>	ID

^{*} Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).
¹ Not all species encountered during clinical performance evaluations

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Lactococcus plantarum</i>	ID	<i>Staphylococcus carnosus</i>	ID/AST
<i>Lactococcus raffinolactis</i>	ID	<i>Staphylococcus chromogenes</i>	ID/AST
<i>Leifsonia aquatica</i>	ID	<i>Staphylococcus cohnii</i>	ID/AST
<i>Leuconostoc citreum</i>	ID	<i>Staphylococcus cohnii ssp cohnii</i>	ID/AST
<i>Leuconostoc lactis</i>	ID	<i>Staphylococcus cohnii ssp urealyticum</i>	ID/AST
<i>Leuconostoc mesenteroides ssp cremoris</i>	ID	<i>Staphylococcus condimentii</i>	AST
<i>Leuconostoc mesenteroides ssp mesenteroides</i>	ID	<i>Staphylococcus delphini</i>	AST
<i>Leuconostoc pseudomesenteroides</i>	ID	<i>Staphylococcus epidermidis</i>	ID/AST
<i>Listeria grayi</i>	ID	<i>Staphylococcus equorum</i>	ID/AST
<i>Listeria innocua</i>	ID	<i>Staphylococcus felis</i>	ID/AST
<i>Listeria ivanovii</i>	ID	<i>Staphylococcus fleurettii</i>	AST
<i>Listeria monocytogenes</i>	ID	<i>Staphylococcus gallinarum</i>	ID/AST
<i>Listeria welshimeri</i>	ID	<i>Staphylococcus haemolyticus</i>	ID/AST
<i>Lysinibacillus sphaericus</i>	ID	<i>Staphylococcus hominis</i>	ID/AST
<i>Macrococcus caseolyticus</i>	ID	<i>Staphylococcus hyicus</i>	ID/AST
<i>Micrococcus luteus</i>	ID	<i>Staphylococcus intermedius</i>	ID/AST
<i>Micrococcus lylae</i>	ID	<i>Staphylococcus kloosii</i>	ID/AST
<i>Oerskovia xanthineolytica</i>	ID	<i>Staphylococcus lentus</i>	ID/AST
<i>Paenibacillus alvei</i>	ID	<i>Staphylococcus lugdunensis</i>	ID/AST
<i>Paenibacillus macerans</i>	ID	<i>Staphylococcus lutrae</i>	AST
<i>Pediococcus acidilactici</i>	ID	<i>Staphylococcus muscae</i>	AST
<i>Pediococcus damnosus</i>	ID	<i>Staphylococcus pasteurii</i>	ID/AST
<i>Pediococcus dextrinicus</i>	ID	<i>Staphylococcus piscifermentans</i>	AST
<i>Pediococcus parvulus</i>	ID	<i>Staphylococcus pulvereri</i>	AST
<i>Pediococcus pentosaceus</i>	ID	<i>Staphylococcus saccharolyticus</i>	AST
<i>Rhodococcus equi</i>	ID	<i>Staphylococcus saprophyticus</i>	ID/AST
<i>Rothia dentocariosa</i>	ID	<i>Staphylococcus schleiferi</i>	ID/AST
<i>Rothia mucilaginosa</i>	ID	<i>Staphylococcus schleiferi ssp coagulans</i>	ID/AST
<i>Staphylococcus arlettae</i>	AST	<i>Staphylococcus schleiferi ssp schleiferi</i>	ID/AST
<i>Staphylococcus aureus</i>	ID/AST	<i>Staphylococcus sciuri</i>	ID/AST
<i>Staphylococcus auricularis</i>	ID/AST	<i>Staphylococcus simulans</i>	ID/AST
<i>Staphylococcus capitis</i>	ID/AST	<i>Staphylococcus succinus ssp casei</i>	AST
<i>Staphylococcus capitis ssp capitis</i>	ID/AST	<i>Staphylococcus succinus ssp succinus</i>	AST
<i>Staphylococcus capitis ssp ureolyticus</i>	ID/AST	<i>Staphylococcus vitulinus</i>	ID/AST
<i>Staphylococcus caprae</i>	ID/AST		

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Staphylococcus warneri</i>	ID/AST	<i>Streptococcus equinus</i>	ID
<i>Staphylococcus xylosus</i>	ID/AST	<i>Streptococcus gordonii</i>	ID
<i>Streptococcus acidominimus</i>	ID	<i>Streptococcus</i> group C/G (large colony)	ID
<i>Streptococcus agalactiae</i> (Strep. group B)	ID	<i>Streptococcus intermedius</i>	ID
<i>Streptococcus anginosus</i>	ID	<i>Streptococcus mitis</i>	ID
<i>Streptococcus bovis</i> (Strep. group D)	ID	<i>Streptococcus mutans</i>	ID
<i>Streptococcus bovis</i> I (Strep. group D)	ID	<i>Streptococcus oralis</i>	ID
<i>Streptococcus bovis</i> II (Strep. group D)	ID	<i>Streptococcus parasanguinis</i>	ID
<i>Streptococcus canis</i>	ID	<i>Streptococcus pneumoniae</i>	ID
<i>Streptococcus constellatus</i>	ID	<i>Streptococcus porcinus</i>	ID
<i>Streptococcus cristatus</i>	ID	<i>Streptococcus pyogenes</i> (Strep. group A)	ID
<i>Streptococcus dysgalactiae</i> ssp <i>dysgalactiae</i>	ID	<i>Streptococcus salivarius</i>	ID
<i>Streptococcus dysgalactiae</i> ssp <i>equisimilis</i>	ID	<i>Streptococcus sanguinis</i>	ID
<i>Streptococcus equi</i>	ID	<i>Streptococcus sobrinus</i>	ID
<i>Streptococcus equi</i> ssp <i>equi</i>	ID	<i>Streptococcus uberis</i>	ID
<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>	ID	<i>Streptococcus vestibularis</i>	ID

F.3.4 Gram Positive (0.25 McFarland)

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Aerococcus urinae</i>	ID	<i>Enterococcus raffinosus</i>	ID/AST
<i>Aerococcus viridans</i>	ID	<i>Gemella haemolysans</i>	ID
<i>Alloiococcus otitidis</i>	ID	<i>Gemella morbillorum</i>	ID
<i>Dermaecoccus nishinomiyaensis</i>	ID	<i>Globicatella sanguinis</i>	ID
<i>Enterococcus avium</i>	ID/AST	<i>Helcococcus kunzii</i>	ID
<i>Enterococcus casseliflavus</i>	ID/AST	<i>Kocuria kristinae</i>	ID
<i>Enterococcus durans</i>	ID/AST	<i>Kocuria rosea</i>	ID
<i>Enterococcus faecalis</i>	ID/AST	<i>Kocuria varians</i>	ID
<i>Enterococcus faecium</i>	ID/AST	<i>Kytococcus sedentarius</i>	ID
<i>Enterococcus gallinarum</i>	ID/AST	<i>Lactococcus lactis</i> ssp <i>cremoris</i>	ID
<i>Enterococcus hirae</i>	ID/AST	<i>Lactococcus lactis</i> ssp <i>hordniae</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).
¹ Not all species encountered during clinical performance evaluations

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Lactococcus plantarum</i>	ID	<i>Staphylococcus lugdunensis</i>	ID/AST
<i>Leuconostoc citreum</i>	ID	<i>Staphylococcus pasteurii</i>	ID/AST
<i>Leuconostoc lactis</i>	ID	<i>Staphylococcus saprophyticus</i>	ID/AST
<i>Leuconostoc mesenteroides</i> ssp <i>mesenteroides</i>	ID	<i>Staphylococcus schleiferi</i> ssp <i>coagulans</i>	ID/AST
<i>Listeria innocua</i>	ID	<i>Staphylococcus schleiferi</i> ssp <i>schleiferi</i>	ID/AST
<i>Listeria monocytogenes</i>	ID	<i>Staphylococcus sciuri</i>	ID/AST
<i>Macrococcus caseolyticus</i>	ID	<i>Staphylococcus simulans</i>	ID/AST
<i>Micrococcus luteus</i>	ID	<i>Staphylococcus vitulinus</i>	ID/AST
<i>Micrococcus lylae</i>	ID	<i>Staphylococcus warneri</i>	ID/AST
<i>Pediococcus acidilactici</i>	ID	<i>Staphylococcus xylosus</i>	ID/AST
<i>Pediococcus damnosus</i>	ID	<i>Streptococcus agalactiae</i> (Strep. group B)	ID
<i>Pediococcus dextrinicus</i>	ID	<i>Streptococcus anginosus</i>	ID
<i>Pediococcus parvulus</i>	ID	<i>Streptococcus bovis</i> (Strep. group D)	ID
<i>Pediococcus pentosaceus</i>	ID	<i>Streptococcus bovis</i> I (Strep. group D)	ID
<i>Rothia mucilaginosa</i>	ID	<i>Streptococcus bovis</i> II (Strep. group D)	ID
<i>Staphylococcus aureus</i>	ID/AST	<i>Streptococcus constellatus</i>	ID
<i>Staphylococcus auricularis</i>	ID/AST	<i>Streptococcus cristatus</i>	ID
<i>Staphylococcus capitis</i>	ID/AST	<i>Streptococcus equi</i>	ID
<i>Staphylococcus caprae</i>	ID/AST	<i>Streptococcus gordonii</i>	ID
<i>Staphylococcus carnosus</i>	ID/AST	<i>Streptococcus</i> group C/G (large colony)	ID
<i>Staphylococcus chromogenes</i>	ID/AST	<i>Streptococcus intermedius</i>	ID
<i>Staphylococcus cohnii</i> ssp <i>cohnii</i>	ID/AST	<i>Streptococcus mitis</i>	ID
<i>Staphylococcus cohnii</i> ssp <i>urealyticum</i>	ID/AST	<i>Streptococcus mutans</i>	ID
<i>Staphylococcus epidermidis</i>	ID/AST	<i>Streptococcus oralis</i>	ID
<i>Staphylococcus equorum</i>	ID/AST	<i>Streptococcus parasanguinis</i>	ID
<i>Staphylococcus felis</i>	ID/AST	<i>Streptococcus pneumoniae</i>	ID
<i>Staphylococcus gallinarum</i>	ID/AST	<i>Streptococcus porcinus</i>	ID
<i>Staphylococcus haemolyticus</i>	ID/AST	<i>Streptococcus pyogenes</i> (Strep. group A)	ID
<i>Staphylococcus hominis</i>	ID/AST	<i>Streptococcus salivarius</i>	ID
<i>Staphylococcus hyicus</i>	ID/AST	<i>Streptococcus sanguinis</i>	ID
<i>Staphylococcus intermedius</i>	ID/AST	<i>Streptococcus sobrinus</i>	ID
<i>Staphylococcus kloosii</i>	ID/AST	<i>Streptococcus uberis</i>	ID
<i>Staphylococcus lentus</i>	ID/AST	<i>Streptococcus vestibularis</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).
¹ Not all species encountered during clinical performance evaluations

F.3.5 Streptococci

STREPTOCOCCI TAXA ¹	ID, AST, ID/AST*	STREPTOCOCCI TAXA ¹	ID, AST, ID/AST*
<i>Streptococcus acidominimus</i>	ID/AST	<i>Streptococcus</i> group C/G (large colony)	ID/AST
<i>Streptococcus agalactiae</i> (Strep. group B)	ID/AST	<i>Streptococcus</i> group CFG (small colony)	AST
<i>Streptococcus alactolyticus</i>	AST	<i>Streptococcus</i> group D (non-enterococcus)	AST
<i>Streptococcus anginosus</i>	ID/AST	<i>Streptococcus</i> group E	AST
<i>Streptococcus anginosus</i> (previously <i>milleri</i>) group	ID/AST	<i>Streptococcus</i> group F	AST
<i>Streptococcus bovis</i> (Strep. group D)	AST	<i>Streptococcus</i> group G (large colony)	AST
<i>Streptococcus bovis</i> I (Strep. group D)	ID/AST	<i>Streptococcus</i> group L	AST
<i>Streptococcus bovis</i> II (Strep. group D)	ID/AST	<i>Streptococcus hyointestinalis</i>	AST
<i>Streptococcus canis</i>	ID/AST	<i>Streptococcus infantarius ssp coli</i>	AST
<i>Streptococcus constellatus</i>	ID/AST	<i>Streptococcus infantarius ssp infantarius</i>	AST
<i>Streptococcus criceti</i>	AST	<i>Streptococcus infantis</i>	AST
<i>Streptococcus cristatus</i>	ID/AST	<i>Streptococcus iniae</i>	AST
<i>Streptococcus downei</i>	AST	<i>Streptococcus intermedius</i>	ID/AST
<i>Streptococcus dysgalactiae</i>	AST		
<i>Streptococcus dysgalactiae ssp dysgalactiae</i>	ID/AST	<i>Streptococcus milleri</i> group	AST
<i>Streptococcus dysgalactiae ssp equisimilis</i>	ID/AST	<i>Streptococcus mitis</i>	ID/AST
<i>Streptococcus dysgalactiae/canis</i>	ID/AST	<i>Streptococcus mitis</i> group	ID/AST
<i>Streptococcus equi</i>	ID/AST	<i>Streptococcus mutans</i>	ID/AST
<i>Streptococcus equi ssp equi</i>	ID/AST	<i>Streptococcus mutans</i> group	AST
<i>Streptococcus equi ssp zooepidemicus</i>	ID/AST	<i>Streptococcus oralis</i>	ID/AST
<i>Streptococcus equinus</i>	ID/AST	<i>Streptococcus parasanguinis</i>	ID/AST
<i>Streptococcus ferus</i>	AST	<i>Streptococcus peroris</i>	AST
<i>Streptococcus gallolyticus ssp macedonicus</i>	AST	<i>Streptococcus pneumoniae</i>	ID/AST
<i>Streptococcus gordonii</i>	ID/AST	<i>Streptococcus porcinus</i>	ID/AST
<i>Streptococcus</i> group A (Strep. <i>pyogenes</i>)	AST	<i>Streptococcus pyogenes</i> (Strep. group A)	ID/AST
<i>Streptococcus</i> group B (Strep. <i>agalactiae</i>)	AST	<i>Streptococcus rattii</i>	AST
<i>Streptococcus</i> group C (large colony)	AST	<i>Streptococcus salivarius</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

STREPTOCOCCI TAXA ¹	ID, AST, ID/AST [*]	STREPTOCOCCI TAXA ¹	ID, AST, ID/AST [*]
<i>Streptococcus salivarius</i> group	AST	<i>Streptococcus thermophilus</i>	AST
<i>Streptococcus sanguinis</i>	ID/AST	<i>Streptococcus uberis</i>	ID/AST
<i>Streptococcus sanguinis</i> group	AST	<i>Streptococcus vestibularis</i>	ID/AST
<i>Streptococcus sobrinus</i>	ID/AST	<i>Streptococcus viridans</i> group	AST
<i>Streptococcus suis</i>	AST		

^{*} Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 2.3.1).

¹ Not all species encountered during clinical performance evaluations

F.3.5 Yeast

YEAST TAXA	ID, AST, ID/AST	YEAST TAXA	ID, AST, ID/AST
<i>Blastoschizomyces capitatus</i>	ID	<i>Candida rugosa</i>	ID
<i>Candida albicans</i>	ID	<i>Candida sake</i>	ID
<i>Candida apicola</i>	ID	<i>Candida sphaerica</i>	ID ¹
<i>Candida boidinii</i>	ID ²	<i>Candida tropicalis</i>	ID
<i>Candida bracarensis</i>	ID ¹	<i>Candida utilis</i>	ID
<i>Candida catenulata</i>	ID	<i>Candida viswanathii</i>	ID
<i>Candida ciferrii</i>	ID	<i>Candida zeylanoides</i>	ID
<i>Candida dubliniensis</i>	ID	<i>Cryptococcus albidus</i>	ID
<i>Candida firmetaria</i>	ID	<i>Cryptococcus humicolus</i>	ID
<i>Candida freyschussii</i>	ID	<i>Cryptococcus laurentii</i>	ID ²
<i>Candida glabrata</i>	ID	<i>Cryptococcus luteolus</i>	ID
<i>Candida guilliermondii</i>	ID	<i>Cryptococcus neoformans</i>	ID
<i>Candida guilliermondii</i> var <i>membranaefaciens</i>	ID	<i>Cryptococcus terreus</i>	ID ¹
<i>Candida haemulonii</i>	ID	<i>Cryptococcus uniguttulatus</i>	ID
<i>Candida inconspicua</i>	ID	<i>Exophiala</i> species	ID
<i>Candida kefyr</i>	ID	<i>Geotrichum</i> species	ID
<i>Candida krusei</i>	ID	<i>Hortaea werneckii</i>	ID
<i>Candida lipolytica</i>	ID	<i>Kloeckera</i> species	ID ¹
<i>Candida lusitanae</i>	ID	<i>Malassezia furfur</i> complex	ID ¹
<i>Candida magnoliae</i>	ID	<i>Malassezia pachydermatis</i>	ID ¹
<i>Candida melibiosica</i>	ID	<i>Malassezia sympodialis</i>	ID ¹
<i>Candida membranaefaciens</i>	ID	<i>Pichia burtonii</i>	ID
<i>Candida norvegensis</i>	ID	<i>Pichia farinosa</i>	ID
<i>Candida parapsilosis</i> complex	ID	<i>Prototheca wickerhamii</i>	ID
<i>Candida pararugosa</i>	ID	<i>Prototheca zopfii</i>	ID
<i>Candida pelliculosa</i>	ID	<i>Rhodotorula glutinis</i>	ID
<i>Candida pulcherrima</i>	ID	<i>Rhodotorula minuta</i>	ID

¹ With SABDX, SABEM, SABHI media only

² With CHOC, COLSB, TSASB media only

YEAST TAXA	ID, AST, ID/AST	YEAST TAXA	ID, AST, ID/AST
<i>Rhodotorula mucilaginosa</i> var <i>mucilaginosa</i>	ID	<i>Trichosporon loubieri</i>	ID ¹
<i>Saccharomyces cerevisiae</i>	ID	<i>Trichosporon mucoides</i>	ID
<i>Sporobolomyces salmonicolor</i>	ID	<i>Trichosporon ovoides</i>	ID ²
<i>Trichosporon asahii</i>	ID	<i>Wangiella dermatitidis</i>	ID
<i>Trichosporon inkin</i>	ID ¹	<i>Zygosaccharomyces bailli</i>	ID ¹

¹ With SABDX, SABEM, SABHI media only

¹ With CHOC, COLSB, TSASB media only

Bibliography

- 1 Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. *Am. J. Public. Health.* 8:922-923.
- 2 Arnold, W.M., Jr., and R.H. Weaver. 1948. Quick microtechniques for the identification of cultures. I. Indole production. *J. Lab. Clin. Med.* 33:195-195.
- 3 Bachmann, B., and R.H. Weaver. 1951. Rapid microtechniques for identification of cultures. V. Reduction of nitrates. *Am. J. Clin. Pathol.* 21:195-196.
- 4 Hannan, J., and R.H. Weaver. 1949. Quick microtechniques for the identification of cultures II. Germentations. *J. Lab. Clin. Med.* 33:1338-1341.
- 5 Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
- 6 Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. *Appl. Microbiol.* 5:36-40.
- 7 Snyder, M.L. 1954. Paper discs containing entire culture medium for the differentiation of bacteria. *Pathol. Bacteriol.* 67:217-226.
- 8 Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. *Puerto Rican J. Publ. Hlth. Trop. Med.* 25:96-100.
- 9 Weaver, R.H. 1954. Quicker bacteriological results. *Am. J. Med. Technol.* 20:14-26.
- 10 Kampfer, P., O. Rauhoff, And W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. *J. Clin. Microbiol.* 29:2877-2879.
- 11 Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. *Microbiol. Rev.* 55:335-348.
- 12 Rammelkamp, C.H., and T. Maxon. 1942. Resistance of *Staphylococcus aureus* to the action of penicillin. *Proc. Soc. Biol. And Med.* 51:386-389.
- 13 Marymont, J.H. and R.M. Wentz. 1966. Serial dilution antibiotic sensitivity testing with the microtiter system. *Am. J. Clin. Pathol.* 45:548-551.

- 14 Gavan, R.L., and M.A. Town. 1970. A microdilution method for antibiotic susceptibility testing: an evaluation. *Am. J. Clin. Pathol.* 53:880-885.
- 15 Lancaster, M.V. and D. Rebecca. 1996. Antibiotic and cytotoxic drug susceptibility assays using resazurin and poisoning agents. US Patent #5,501,959.
- 16 Data on file in Becton Dickinson Diagnostics.
- 17 Murray, Patrick R., et al. ed., *Manual of Clinical Microbiology*, 7th Edition, ASM Press, Washington, D.C., 1999.
- 18 DM Livermore. 1995. Beta-Lactamases in laboratory and clinical resistance. *Clin. Microbiol. Rev.* 8: 557-584.
- 19 K Bush, GA Jacoby, and AA Medeiros. 1995. A functional classification scheme for beta-lactamases and its correlation with molecular structure. *Antimicrob. Agents Chemother.* 39: 1211-1233.
- 20 GA Jacoby and P Han. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. *J. Clin. Microbiol.* 34: 908-911.
- 21 Kenneth S. Thomson, Christine C. Sanders, and Ellen Smith Moland. 1999. Use of Microdilution Panels with and without -Lactamase Inhibitors as a Phenotypic Test for -Lactamase Production among *Escherichia coli*, *Klebsiella* spp., *Enterobacter* spp., *Citrobacter freundii*, and *Serratia marcescens*. *Antimicrob. Agents Chemother.* 43: 1393-1400.
- 22 CLSI. M100-S18 *Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement*. January, 2008.
- 23 CLSI. M07-A8 *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard-Eighth Edition*. January, 2009.
- 24 Livermore, D.M., T.G. Winstanley, and K.P. Shannon. 2001. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. *J. Antimicrob. Chemother.* 48, Suppl. S1, 87-102.
- 25 Cauwelier, et al. "Evaluation of a disk diffusion method with cefoxitin (30 mcg) for detection of methicillin-resistant *Staphylococcus aureus*." *Eur J. Clin. Microbiol Infect Dis.* May 2004,23(5):389-392.
- 26 Barnett, J.A., Payne, R.W. and Yarrow, D., *Yeasts: Characteristics and Identification, 2nd Edition*, Cambridge University Press, Cambridge, U.K. 1990.
- 27 Larone, D. H., *Medically Important Fungi, A Guide to Identification, 4th Edition*, ASM Press, Washington, D.C. 2002
28. de Hoog, G.S., Guarro, J., Gené, J. and Figueras, M.J., *Atlas of Clinical Fungi, 2nd Edition*, Centraalbureau voor Schimmelcultures, Utrecht, the Netherlands/Universita Rovira I Virgili, Reus, Spain, 2000.

Event Log Messages

J.1 LIS Related Messages

The View Syslog function, accessible through the Maintenance Menu (see Section 6.2.4.13), enables you to see messages that the instrument writes to the Event Log. The messages are in the following format:

date time message type: message text

where *date* represents the day, month, and year in the format you have chosen

time represents the time in the format you have chosen

message type is one of the following:

- LIS Interface Message
- LIS Unsolicited Message
- LIS Order Cancelled
- LIS Query Assembly
- LIS Receiving Query
- LIS Configuration Change

message text is the actual text message that appears

Below is a list of the messages, along with a description of the message and any actions you may perform to correct the problem. These messages are grouped by the *message type* listed above.

LIS Interface Messages

LIS Interface Messages are library messages generated by the LIS manager. They are listed below in alphabetical order.

Bad Frame Received From LIS

DESCRIPTION – This error is generated when the LIS downloads a frame that has not been properly formatted. The error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – The solution to this problem is to review the information that is being sent to the instrument via a communications line monitor. Compare the information captured with the specifications found in the BD – LIS Vendor Interface Document. Any discrepancies observed should be corrected and then the transmission should be attempted again.

Detailed Download Message Attached To Notification

DESCRIPTION – These notifications are generated whenever a complete message is sent or received across the port. The error types are LIS_NOTIFY. The detailed description is not displayed.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Detailed Upload Message Attached To Notification

DESCRIPTION – These notifications are generated whenever a complete message is sent or received across the port. This provides the Host Application with the complete ASTM message string that was exchanged. The error types are LIS_NOTIFY and the detailed descriptions contain the ASTM message string.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Disallowed Characters Contained In Field

DESCRIPTION – This error is generated when an upload field contains characters that are not allowed in an upload message. These characters differ based on logical protocol and will be stripped from the upload field. This error type is LIS_NON_FATAL. The detailed description is not displayed.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Download Field Was Concatenated

DESCRIPTION – This notification is generated when a download field is larger than the maximum size set by the Host Application. The field will be concatenated to the maximum size. This error type is listed as LIS_LOG. The detailed description is not displayed.

CORRECTIVE ACTION(S) – In this case the content that is being assembled needs to be reviewed. Review the BD – Vendor Interface Specification for limitations on field lengths in the message. Then the LIS code should be updated so that the information sent to the instrument remains within the defined limits.

LIS Interface Messages (cont.)

Download Message Has A Bad Or Missing Header Record

DESCRIPTION – This error is generated when a download message does not have a properly formatted header record. The LIS IM cannot continue processing the message and it will be deleted. This error type is listed as LIS_NON_FATAL, and the detailed description holds the ASTM message received from the LIS.

CORRECTIVE ACTION(S) – Check the header line for errors. Correct the header error and resend the message.

Download Physical Communication With LIS Has Begun

DESCRIPTION – This notification is generated when the download thread in the physical layer is actively downloading a message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Download Physical Communication With LIS Has Completed

DESCRIPTION – This notification is generated when the download thread in the logical layer has completed processing a download message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Download Record Is Out Of Sequence Or Has Bad Sequence Number

DESCRIPTION – This error is generated when a download record contains records with a bad sequence number, or violates the record hierarchy for the logical level protocol. This will cause the message to be deleted from the download queue. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Appropriate updates should be made to the application generating the logical message so the structure conforms to the ASTM specifications for the interface. Most questions concerning the message content can be addressed by referencing the BD – Vendor Interface Specifications and the ASTM specifications for the interface.

Early Termination Of Transfer Session By LIS

DESCRIPTION – This error is generated when there is an error in LIS communication. The control character to end a transfer session was received before the appropriate number of characters for the frame were received.

CORRECTIVE ACTION(S) – Review the frames being exchanged between the instrument and the LIS to assure that the proper number of characters are contained within the packets that are being exchanged. If the correct number of characters is not in the packet, update the code so that the correct number of characters is included in each frame.

LIS Interface Messages (cont.)

Expected Frame Not Sent

DESCRIPTION – This error is generated when there is an error in LIS communication. The LIS initiated a transfer session but did not send any data before a time-out occurred. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Connect a communications line monitor and restart the transmission. Observe the transmission to determine if the complete transmission occurs or if the transmission is interrupted early. If the transmission appears to be complete or the LIS appears to be attempting to send the transmission contact your local BD representative for assistance. If the transmission appears to be incomplete from the LIS side of the transmission then the LIS code should be investigated for problems that could terminate the transmission early.

Internal Assert Condition Found In LIS Library

DESCRIPTION – This message indicates that the Instrument has encountered an error that should not occur. This error will result in the instrument rebooting.

CORRECTIVE ACTION(S) – Contact your local BD representative. At the time the problem occurs the application writes data to the log file that indicates the nature of the Assert. The BD representative should instruct you on the appropriate procedure to collect the data for the condition.

LIS Debug Error

DESCRIPTION – This error is generated when there has been a problem internally to the instrument physical layer. This error may be triggered for an index out of bounds, or unsupported memory area IDs, etc. This error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – This error is a Fatal Error. If this error appears in the log file, contact your local BD representative. The sequence of events that produced this message should be documented so that BD representatives can reproduce the error and then provide an appropriate course of action to address the problem.

LIS Interface Messages (cont.)

LIS Did Not Acknowledge Sent Frame

DESCRIPTION – This error is generated when the instrument has sent a frame to the LIS but has not received an acknowledgment before a time-out occurred. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – This message should not occur during normal operation of the instrument interface. If this error is being encountered during development of the interface then the LIS development group should connect a data communication monitor to the serial interface cable and review the information that is being exchanged between the two devices. It is likely that the LIS is not generating the appropriate response to the message that the instrument has sent.

LIS Is Not Responding To Output Request

DESCRIPTION – This error is generated when the instrument is trying to establish a transfer session but the LIS is not responding. When this error is sent, the instrument is assuming the LIS connection is broken or the LIS is down. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Check to assure that the LIS interface is active and ready to receive messages from the instrument. If the LIS is operating correctly, connect a communications line monitor and review the information that has been captured. If the instrument is attempting to establish communications with the LIS, it would be appropriate to review the code that has been written to interface with the instrument and assure that the code is appropriate to receive the information sent by the instrument.

LIS Never Completed Current Frame

DESCRIPTION – This error is generated when there is an error in LIS communication. The control characters expected to end a transmitted frame were never received.

CORRECTIVE ACTION(S) – Review the frames being exchanged between the instrument and the LIS via line monitor. Be sure to review all frames included in the transmission. If any frames do not contain the appropriate termination characters make appropriate changes to the interface code to correct the problem.

Logical Processing Of LIS Data Has Begun

DESCRIPTION/CORRECTIVE ACTION(S) – Message is informational – no action required.

Logical Processing Of LIS Data Has Completed

DESCRIPTION/CORRECTIVE ACTION(S) – Message is informational – no action required.

Message Packet Passed To Host Application

DESCRIPTION – This notification is generated when the LIS IM has successfully passed a message on to the Host Application. This notification exists mostly as a debug message. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – No action required by user. Information is included to help LIS manufacturers debug and implement their interface.

LIS Interface Messages (cont.)

Message Received By LIS

DESCRIPTION – This notification is generated after the LIS has received and acknowledged a message from the instrument.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Message Received From LIS

DESCRIPTION – This notification is generated when a download message has been properly received by the Physical Interface from the LIS. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – This message is an informational message and requires no action to be taken. The transmission was successful. This message will be generated even if the message contained content errors.

Must Re-send Output Frame To LIS

DESCRIPTION – This notification is generated when a frame sent by the LIS IM was not properly received by the LIS. The message frame will be resent, according to low level protocol specifications. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – If this message appears infrequently and the messages that are being exchanged are completed correctly, it is likely that no action is required. If the message above is encountered frequently, there is likely a problem in the interface. The LIS manufacturer and BD representative should be contacted to diagnose the interface to assure that it is operating properly. When the link is operating optimally this message should not appear in the log.

No Queue Memory For Download Messages

DESCRIPTION – The Phoenix Instrument is designed to operate in a limited amount of memory (as defined in the configuration structure). When these memory resources are full with download or upload messages, this notification will be generated. The notification may be common if the LIS tries to download too much information at once. The Phoenix Instrument will NAK data for which it cannot allocate memory, and the LIS will have to resend the data.

CORRECTIVE ACTION(S) – Resend the messages that were rejected by the instrument.

No Queue Memory For Upload Messages

DESCRIPTION – The Phoenix Instrument is designed to operate in a limited amount of memory (as defined in the configuration structure). When these memory resources are full with download or upload messages, this notification will be generated. This message should not be encountered when using the instrument LIS interface.

CORRECTIVE ACTION(S) – If this message is encountered, the problem should be reported to the BD representative. It will be helpful if the process that generated the problem is documented.

LIS Interface Messages (cont.)

No Response Received From Previous LIS Request Message

DESCRIPTION – This notification is generated when the LIS IM has uploaded a query but no response is received. This notification indicates that the original query was cancelled. This notification should occur some time after the query, as determined by the value in the configuration structure. The error type is listed as LIS_LOG. The detailed description is not displayed in the current version.

CORRECTIVE ACTION(S) – Validate that the LIS system is correctly connected to the instrument and that the LIS interface has been activated.

Operating System Error

DESCRIPTION – This error is generated when the operating system class encounters an error in one of its routines. This could be caused by a number of OS errors, including not properly initializing a port or not properly creating an event handle, and others. These errors should not occur under normal conditions and the LIS IM cannot recover from them. The error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Output Message Was Sent To LIS

DESCRIPTION – This notification is generated when an upload message has been successfully transmitted to the LIS. This error type is listed as LIS_NOTIFY. Notification sent when the instrument begins sending a message to the LIS. Receiving this notification does not signify the message was accepted by the LIS.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Queue Memory For Download Messages Free

DESCRIPTION – This notification is used in pairs with the No Queue Memory For Download Messages above. When memory has previously been determined to be full, and now has been released, this notification will be generated.

CORRECTIVE ACTION(S) – This is an indication to a user of a batch-oriented interface that the instrument is capable of receiving the next group of messages.

Queue Memory For Upload Messages Free

DESCRIPTION – This notification is used in pairs with the No Queue Memory For Upload Messages above. When memory has previously been determined to be full, and now has been released, this notification will be generated.

CORRECTIVE ACTION(S) – This message should not appear in the log during operation of the interface. If this message appears in the log, contact your local BD representative.

LIS Interface Messages (cont.)

Response Message Received

DESCRIPTION – This notification is generated when the LIS has downloaded a response to a BD generated query. This is also an indication that the LIS IM can upload another query to the LIS. The error type is LIS_NOTIFY. The detailed description is not displayed in the current version.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Unsupported Field In Configuration File

DESCRIPTION – This error is generated during startup when a field in the configuration structure does not match one of the supported configurations of the LIS IM. This error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Upload Physical Communication With LIS Has Begun

DESCRIPTION – This notification is generated when the upload thread in the physical layer is actively uploading a message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Upload Physical Communication With LIS Has Completed

DESCRIPTION – This notification is generated when the upload thread in the physical layer has completed sending a message to the LIS.

CORRECTIVE ACTION(S) – Message is informational – no action required.

LIS Unsolicited Message

The following sections contain the possible messages associated with Event Log entries in response to unsolicited download requests. They are listed below in alphabetical order.

LIS Response Error

DESCRIPTION – This message is set after the LIS sends an incorrectly formatted query message to the instrument. We will respond by sending a response error message to the LIS.

CORRECTIVE ACTION(S) – If the interface is operating correctly, this message should not be present in the log. If this message is found in the log during validation or development of the interface, the code responsible for generating the queries should be reviewed and adjusted so that the query is assembled without the invalid field. The invalid field will be the last "LIS Query Assembly" message listed before this message in the log file.

LIS Sent Query Before Previous Query Completed

DESCRIPTION – This message is set if the LIS initiates a query and while the instrument is busy with this query, the LIS cancels the first query and immediately starts another.

CORRECTIVE ACTION(S) – If the interface is operating correctly this message should not be present in the log. If this message is found in the log during validation or development of the interface, the code responsible for generating the queries should be reviewed and adjusted so that it will not send another query while it is waiting for an outstanding query to be canceled.

LIS Order Cancelled Messages

(Message consists of a field name shown below at * and invalid field contents.)

DESCRIPTION – The validation for the noted field failed.

- * Invalid Sequence Number Field
- Missing Accession Number Field
- Invalid Accession Number Field [cannot be 12 digits and begin with 42 or 50 – 59; see Section 5 – Reference for other requirements]
- Missing Isolate Field
- Invalid Isolate Field
- Invalid Panel Usage
- Invalid Test ID Field
- Invalid Test Strain for this QC Panel Type
- Mismatch Sequence Number & Test ID
- Missing Test ID Field
- Invalid Organism ID Field
- Invalid Priority Field
- Problem Storing Record to Database

When an error is detected, the message for that error is logged and the remainder of the message checking is terminated. There could be other errors in the order that were not reported. The error checking priority depends upon whether the order contains a sequence number or not.

CORRECTIVE ACTION(S) – In an operational interface it is possible to encounter these errors if the LIS implementation allows entry of data into the fields that violate the Instrument field rules. In this case the user of the LIS should limit the characters entered into the fields sent to the instrument to those characters noted in the System User's Manual as valid for the field. Once the information has been updated to conform with the instrument's rules the order should be sent to the instrument again.

If the interface is under development the developer of the interface should limit the data being sent to the instrument to those characters that are appropriate for the fields.

LIS Query Assembly Messages

All of the following messages are written to the log file regardless of whether they are invalid or valid. They are placed in the log so that they can be used in conjunction with the messages that follow them. If an error is found in one of the fields, the contents of the field in error will be logged.

(Message consists of a field name shown below at **.)

DESCRIPTION – The message is logged when the instrument uploads a query to the LIS.

** Accession Number
 Sequence Number
 Test ID
 Test Status
 Result Qualifier
 Time Qualifier
 Start Date/Time
 End Date/Time

CORRECTIVE ACTION(S) – If the query is valid, the message is informational – no action required. If the query is invalid, then review the field in error, correct, and resend the query.

(Message consists of a field name and value shown below at ***)

DESCRIPTION – This message is logged when the instrument uploads a query to the LIS.

*** Sequence Number and the value

CORRECTIVE ACTION(S) – Message is informational – no action required.

LIS Configuration Change Messages

(Message consists of a configuration value shown below at ** with old and new values.)

DESCRIPTION – The message is logged when the LIS configuration is changed.

- ** LIS Enabled Value Changed
- Send Interpretation Results Option Changed
- Unsolicited Queries Option Changed
- Send When Placed In Instrument Option Changed
- Results Upload Options Changed
- Baud Value Changed
- Data Bits Value Changed
- Parity Value Changed
- Stop Bits Value Changed
- Packed Frames Value Changed

CORRECTIVE ACTION(S) – Message is informational – no action required.

Notes

Notes section with horizontal lines for writing.

Notes

[illegible]

Index

A

access levels 2-13, 2-14
alarms 4-45
antimicrobial
 used in Phoenix panels F-1, F-2, F-3, F-4, F-5,
 F-6
AST Broth (also AST-S Broth)
 storage and handling 4-2
 use of 4-7, 4-8, 4-14, 4-15
AST Indicator Solution (also AST-S Indicator
Solution)
 storage and handling 4-2
 use of 4-8, 4-9, 4-10, 4-14, 4-15
AST results 4-35, 4-39, 5-17, 5-18, 5-19
 Resistance Markers 2-37, 3-8, 4-40
audible alarms
 meaning 3-8
 setting the volume 2-8
 verification 6-4

B

barcode scanners 3-6
 cleaning 6-7
 replacement 6-46
BD Diagnostics
 international contacts D-1

BD Phoenix system
 alarms and errors 4-45
 cautions 1-14, 1-15
 daily maintenance 6-2
 decontamination 6-8
 explanation of the test 1-2
 installation 2-2
 instrument layout 3-2
 instrument overview 1-6, 1-7, 1-8
 intended use 1-1
 panel overview 1-9
 panel testing 4-26
 precautions 1-6
 principles of the procedure 1-2, 1-3, 1-4, 1-5
 quality control 4-15, 4-16, 4-17, 4-18, 4-19
 replacement parts B-1
 service 7-1
 software installation 2-43
 software overview 1-8
 software setup 2-4, 2-5, 2-6, 2-7, 2-8, 2-9,
 2-10, 2-11, 2-12, 2-13, 2-14
 software update log C-1
 specifications 2-2, 2-3
 startup process 2-43
 system alerts 5-61
 testing overview 1-10
 warnings 1-14, 1-15
 weekly maintenance 6-4
BDXpert system
 viewing triggered rules 4-40

C

cautions 1-15
critical panels 2-35, 2-36, 3-8, 4-20, 4-22, 4-34

D

date
 setting 2-11
 setting format 2-10
decontamination 6-8
door
 unlocked indicator 3-5

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floppy disk
 eject button 3-6
 indicator 3-6

H

how to
 adjust for daylight savings time 2-10
 configure custom rules 2-26, 2-29, 2-30, 2-31, 7-7
 copy configuration from one instrument to another 6-13
 disable ID or AST test on a panel 4-20, 6-27
 enable/disable BDXpert rules 2-34
 enable/disable LIS communications 2-15
 enter hospital/laboratory name/address 2-12
 inoculate panels 4-8, 4-10
 insert panels in the instrument 4-24
 log in panels 4-19
 log in panels downloaded from a LIS system 4-22
 obtain service 7-1
 prepare panels 4-4, 4-6, 4-8, 4-9, 4-10, 4-11, 4-13
 print reports 4-47
 resolve panels needing attention 4-28
 scan a barcode 3-6
 select rule set 2-8
 set audible alarm volume 2-8
 set date format 2-10
 set inoculum density 2-12, 2-13
 set instrument number 2-7
 set language 2-7
 set LIS codes for antimicrobials 2-22
 set LIS codes for organisms 2-20
 set LIS communications parameters 2-14
 set the date 2-11
 set the time 2-11
 set time format 2-11
 setup the software 2-4, 2-5, 2-6, 2-7, 2-10, 2-11
 test a QC panel 4-15, 4-18
 update panel database 6-23, 6-25
 use the instrument interface 4-19
 view panel inventory 5-58
 view panel results 4-35, 4-37, 4-38

I

icons 2-32, 2-35, 2-39
 BDXpert Triggered Rules screen 4-40, 5-21
 configuration 2-5, 2-6, 2-24, 2-43
 Finalization screen 5-29
 Main Status screen 5-4
 Needs Attention screen 5-60, 5-63
 Panel Inventory screen 5-60
 Panel Login screen 5-6

- Panel Results screen 5-12
- System Alerts screen 5-61
- ID Broth
 - storage and handling 4-2
 - use of 4-6
- ID results 4-32, 4-33, 5-20
 - adding 4-38
 - modifying 4-38
- inoculum density 2-12, 2-13, 4-6, 4-8, 4-11
- instrument
 - automatic light source adjustment 1-11
 - barcode scanners 1-7
 - carousel 1-6
 - control electronics 1-8
 - drive mechanism 1-6
 - incubation 1-7
 - normalizers 1-10
 - optical system 1-6
 - panel status indicators 1-7
 - symbols appearing on 1-13
- interpretation rules
 - activating/deactivating BDXpert system 2-8
 - configuring custom rules 2-26, 2-29, 2-31, 7-7
 - enabling/disabling BDXpert rules 2-34
 - how to select rule set 2-8
 - viewing rule set version 2-8
 - viewing triggered BDXpert rules 4-38

K

- keyboard
 - key functions 3-9
 - location 3-2
- keys
 - DOWN/UP ARROW 3-3
 - LOAD PANELS 3-5
 - SILENCE ALARM 3-5
 - soft(ware) 3-4
 - TAB 3-4
 - UNLOAD PANELS 3-5

L

- LCD display
 - automatic dimming 3-5
- LIS communications
 - operation 4-50, 4-52
 - setup 2-15
 - viewing Syslog 6-31, 6-32, 6-33, 6-34

M

- maintenance
 - air filter replacement 6-4, 6-5, 6-6
 - daily 4-2
 - decontamination 6-8
 - Install Panel Configuration 6-22
 - Modify Panel Usage 6-25
 - Phoenix Update Disk 6-24
 - Read Light Source from Diskette 6-29
 - Restore LIS Codes from Disk 6-14
 - Restore Panel Lot Definitions 6-18
 - Restore User Data from Disk 6-11
 - Save LIS Codes to Disk 6-13
 - Save Panel Configuration 6-20
 - Save Panel Lot Definitions 6-16
 - Save User Data to Disk 6-9
 - software maintenance menu 6-9
 - View Syslog 6-31
 - weekly 6-4
 - who performs 6-1
 - Write Light Source to Disk 6-27
- manual
 - conventions 1-12
 - notes, cautions, and warnings 1-14
 - structure 1-11
 - summary of cautions and warnings 1-15
 - use of 1-12

O

- On/Off switch 3-3

P

- panels
 - antimicrobials F-1, F-2, F-3, F-5, F-6
 - automatic association of ID results 4-33
 - batch finalization 4-42, 4-44
 - critical 2-35, 2-36, 3-8, 4-20, 4-22, 4-34
 - disabling ID or AST test on panel 6-27
 - disposal 4-48
 - finalizing 5-16
 - handling recommendations 4-5, 4-9, 4-11, 4-13
 - how to log in 4-19
 - how to place in instrument 4-24
 - how to test 4-26
 - ID results 4-32, 4-33
 - inoculating 4-8, 4-10
 - inoculation station (photo) 1-9
 - inspection after inoculation 4-8, 4-10, 4-13

- isolates 4-2
- labels (applying) 4-6
- materials required for preparation 4-4, 4-5
- media to avoid 4-4
- microbiological hazards 4-4
- modifying panel usage 6-25
- Needs Attention 4-28
- overview 1-9
- Phoenix Update Disk 6-24
- preparing 4-4, 4-7, 4-8, 4-9, 4-10, 4-11, 4-13
- reagents F-8
- storage and handling 4-2
- taxa for ID/AST determination F-15
- timing requirements 4-5, 4-8, 4-9, 4-10, 4-11, 4-13
- transportation caddy (photo) 1-10
- unloading 4-48
- updating panel database 4-49, 6-25
- viewing results 4-35, 4-37, 4-38
- password administration 2-13, 2-14
- power switch 3-3

Q

- quality control 4-18

R

- rapid reporting of results 2-34, 2-35, 2-36, 2-37
- reagents
 - used in Phoenix panels F-8
- replacement parts B-1
- reports 4-47
 - abbreviated Lab Report 2-9
 - Accession Lab 5-36
 - Antimicrobial Code 5-53
 - BDXpert Rule Set Database 5-51
 - Completed Lab 5-33
 - Cumulative QC 5-43
 - Custom Breakpoint Difference 2-25
 - Daily Instrument 5-45
 - enabling/disabling auto printing 2-11
 - entering hospital/laboratory name/address 2-12
 - Finalization Summary 4-43, 4-44
 - Interpretation Rule Set 5-47
 - Lab Report 5-27
 - Needs Attention List 5-38
 - Organism ID Code List 5-52
 - Panel Lot 5-55
 - Panel Lot Database 5-57
 - QC Lot Number 5-41

- Resident Panel 5-39
- Resistance Markers 2-37, 3-8, 4-35
 - setting auto report time 2-11
- Resistance Markers 2-37, 3-8, 4-36, 4-40

S

- soft keys
 - location 3-4
- software
 - how to use 4-19
 - installation 2-43
 - menu tree 5-2
 - overview 1-8
 - setup 2-4, 2-5, 2-6, 2-7, 2-10, 2-11
 - update log C-1
- stations
 - status indicators 3-7
 - verifying LEDs 6-4
- system alerts
 - "E" codes 7-2
 - "W" codes 7-2
 - indicator 3-5
 - reviewing 7-2
 - sub-codes 7-3
 - types 7-2
 - verification 6-4

T

- time
 - setting 2-11
 - setting daylight savings time 2-10
 - setting format 2-11
 - time zone offset 2-10

W

- W100 7-16
- W101 7-17
- W102 7-17
- W200 7-17
- W201 7-17
- W202 7-17
- W203 7-18
- W300 7-18
- W301 7-18
- W302 7-18
- W303 7-19
- W304 7-19
- W305 7-19
- W401 7-19
- W500 7-20

W501 7-20
W502 7-20
W503 7-20
W504 7-21
W505 7-21
W506 7-21
W507 7-21
W510 7-21
W600 7-22
W601 7-22
W602 7-22
W603 7-22
W604 7-23
W605 7-23
W606 7-23
W607 7-23
W608 7-24
W609 7-24
W700 7-24
W701 7-24
W702 7-24
W800 7-25
W801 7-25
W802 7-25
W803 7-25
warnings 1-15
warranty A-1

Reader Comment Card

BD Phoenix System User's Manual (L-003342)

After you have become familiar with your instrument and its documentation, please take a moment to evaluate the manual. Your comments are a valuable tool for us to improve the usefulness of the manuals.
NOTE: For any evaluation of "average" or below, please add a specific comment on the blank lines at the end.

Thank you!

	Excellent	Good	Average	Fair	Poor
Chapters					
Table of Contents	☐	☐	☐	☐	☐
Introduction	☐	☐	☐	☐	☐
Installation	☐	☐	☐	☐	☐
Controls and Indicators	☐	☐	☐	☐	☐
Operation	☐	☐	☐	☐	☐
Reference	☐	☐	☐	☐	☐
Maintenance	☐	☐	☐	☐	☐
Troubleshooting	☐	☐	☐	☐	☐
Performance and Limitations	☐	☐	☐	☐	☐
Index	☐	☐	☐	☐	☐
Appendices	☐	☐	☐	☐	☐
General					
Organization/order of information	☐	☐	☐	☐	☐
Ease of finding information	☐	☐	☐	☐	☐
Clarity of information	☐	☐	☐	☐	☐
Completeness of information	☐	☐	☐	☐	☐
Accuracy of information	☐	☐	☐	☐	☐
Diagrams – Quantity	☐	☐	☐	☐	☐
Diagrams – Quality	☐	☐	☐	☐	☐
Diagrams – Usefulness	☐	☐	☐	☐	☐
Overall appearance	☐	☐	☐	☐	☐
Explanation of hardware	☐	☐	☐	☐	☐
Explanation of software	☐	☐	☐	☐	☐

Comments _____

Optional Name _____

 Hospital _____

 City, State _____

Fold Under & Tape

Fold Here

BUSINESS REPLY MAIL

PERMIT NO. 25

SPARKS, MD

POSTAGE WILL BE PAID BY ADDRESSEE

BD


Technical Writing, Mail Stop 634
7 Loveton Circle / P.O. Box 999
Sparks, MD 21152.9977

NO POSTAGE
NECESSARY IF
MAILED IN THE
UNITED STATES

BD Phoenix™ M50

Automated Microbiology System User's Manual

 **Becton, Dickinson and Company**
7 Loveton Circle
Sparks, MD 21152 USA

Benex Limited
Pottery Road, Dun Laoghaire
Co. Dublin, Ireland

Australian Sponsor:
Becton Dickinson Pty Ltd.
4 Research Park Drive
Macquarie University Research Park
North Ryde, NSW 2113 Australia

Document Number: 500008940(02)

2016-08

REF 443575

English

Not for sale in USA

Change History

Revision	Sections	Reason
01	All	New product release
02	Front Cover and Section 2.5 Icon Chart	Minor typographical corrections

© 2016 BD. BD, the BD Logo, and all other trademarks are property of Becton, Dickinson and Company.

All rights reserved. No part of this publication may be reproduced, transmitted, transcribed, stored in retrieval systems, or translated into any language or computer language, in any form or by any means, electronic, mechanical, magnetic, optical, chemical, manual, or otherwise, without the prior written permission of BD, 7 Loveton Circle, Sparks, MD, 21152, United States of America

Table of Contents

1 – Introduction

1.1 Intended Use	1-1
1.2 Summary and Explanation of the Test	1-1
1.3 Principles of the Procedure	1-1
1.4 Precautions	1-5
1.5 System Overview	1-5
1.5.1 Instrument Overview	1-5
1.5.1.1 Instrument Optical System and Drive Mechanism	1-5
1.5.1.2 Carousel Assembly	1-6
1.5.1.3 Incubation System	1-6
1.5.1.4 Panel Status and Internal Barcode Scanner Assembly	1-6
1.5.1.5 External Barcode Scanner	1-6
1.5.2 Software and Operation Overview	1-6
1.5.3 Panel Overview	1-8
1.5.4 Testing Overview	1-9
1.5.5 Normalizers	1-10
1.6 Limited Warranty	1-11
1.7 Conventions	1-11
1.7.1 Notes, Cautions, and Warnings	1-11
1.7.2 Summary of Cautions and Warnings	1-12
1.8 Bibliography	1-12

2 – Controls, Icons, Indicators, and Installation

2.1 General	2-15
2.2 System Startup	2-15
2.2.1 Instrument Installation	2-15
2.3 Instrument Specifications	2-17
2.3.1 Specification Chart	2-17
2.3.2 Controls and Indicators	2-18
2.3.2.1 Power Switch and Inputs / Outputs	2-18
2.3.2.2 Connector Plate	2-19
2.3.2.3 All-in-one (AIO) PC controls	2-20
2.3.2.4 Power Controls	2-20
2.3.2.5 Other Indicators	2-21
2.3.2.6 Instrument Alert Indicator	2-21

Table of Contents

2.3.2.7 System Alerts	2-22
2.3.2.8 Audible Tones and Alarms	2-22
2.3.2.9 Instrument Door/Panel Button Status	2-23
2.3.2.10 Station Indicators	2-24
2.3.3 Printer	2-25
2.3.4 Barcode Scanners	2-25
2.3.5 Touchscreen Keyboard	2-25
2.3.5.1 Enabling/Disabling Touchscreen Fields	2-25
2.4 User Access and Management Summary	2-25
2.5 Software Icons	2-27
2.6 Software Setup	2-29
2.7 Isolation Mode	2-29
2.8 LIS Operations	2-31
2.8.1 General	2-31
2.8.2 Results Upload Records	2-32
2.8.3 Results Record	2-32
2.8.4 Important Concepts	2-33
2.8.5 Routine System Operation	2-33

3 – User Interface Operations

3.1 Using the Instrument Interface	3-35
3.2 Status Tab	3-35
3.3 Panel Login Tab	3-35
3.3.1 Logging in Panels	3-37
3.3.2 Inserting Panels in the Instrument	3-38
3.3.3 Unloading and Discarding Panels	3-40
3.4 Results Tab	3-41
3.4.1 Adding/Modifying ID Results	3-47
3.4.2 Modifying AST Results	3-47
3.4.3 Modify Panel Usage Sub-Tab	3-48
3.4.4 AST Results Sub-Tab	3-49
3.4.5 BDxpert Rules Sub-Tab	3-52
3.4.6 Lot Information Sub-Tab	3-53
3.4.7 Instrument Organism ID Sub-Tab	3-54
3.4.8 Special Messages Sub-Tab	3-54
3.4.9 Needs Attention Sub-Tab	3-55

Table of Contents

3.4.10 Lab Report	3-55
3.5 Finalization Tab	3-56
3.6 Needs Attention Tab	3-58
3.7 Inventory Tab	3-61
3.8 Reports Tab	3-62
3.8.1 How to Print Reports	3-63
3.8.2 Completed Lab Report	3-63
3.8.3 Accession Lab Report	3-64
3.8.4 Needs Attention List Report	3-64
3.8.5 Resident Panel Report	3-65
3.8.6 QC Lab Report	3-66
3.8.7 Cumulative QC Report	3-66
3.8.8 Daily Instrument Report	3-66
3.8.9 Interpretation Rule Set Report	3-67
3.8.10 BDxpert Rule Set Database Report	3-68
3.8.11 Organism ID Code List Report	3-68
3.8.12 Antimicrobial Code Report	3-68
3.8.13 Lab Report/QC Lab Report	3-68
3.8.14 Finalization Summary Report	3-68
3.8.15 Custom Breakpoint Difference Report	3-69
3.8.16 Current QC Panel Lot Report	3-69
3.8.17 Historical QC Panel Lot Report	3-69
3.8.18 Panel Lot Report	3-69
3.8.19 Panel Lot Database Report	3-69
3.9 Maintenance Tab	3-69
3.9.1 Maintenance Hardware Functions	3-70
3.9.2 Maintenance Software Functions	3-70
3.9.3 Maintenance LIS Functions	3-75
3.9.4 Maintenance Panel Functions	3-75
3.9.5 Maintenance Field Service Functions	3-76
3.9.6 Maintenance Internal BD Use Functions	3-76
3.9.7 Maintenance Event Log Tab	3-76
3.10 Configuration Tab	3-77
3.10.1 Users Sub-Tab	3-77
3.10.2 System Configuration Sub-Tab	3-77

Table of Contents

3.10.3 Communications Sub-Tab.....	3-81
3.10.4 Organism Configuration Sub-Tab	3-83
3.10.5 Antimicrobial Configuration Sub-Tab	3-83
3.10.6 Custom Interpretation Rule Set Sub-Tab	3-84
3.10.7 BDxpert Rules Sub-Tab.....	3-88
3.10.8 Rapid Reporting Sub-Tab	3-89
3.10.9 Panel Lot Definition Sub-Tab	3-91
3.10.10 Instrument Sub-Tab	3-93

4 – Routine Panel Operation

4.1 General	4-95
4.2 Storage and Handling	4-95
4.3 Preparing Panels.....	4-95
4.3.1 Quality Control	4-108
4.4 Automatic Association of Panels.....	4-110
4.4.1 Overview	4-110
4.4.2 Panel Types	4-110
4.4.3 When Auto Association Is Not Performed.....	4-110
4.4.4 Organism ID Conflict Checking	4-110
4.4.5 Auto Association and Related Actions	4-111
4.4.6 Typical Auto Association Examples	4-111
4.4.7 BDxpert Triggered Rules Screen.....	4-113
4.4.8 Resistance Markers	4-114

5 – Maintenance

5.1 General	5-115
5.2 Routine Maintenance	5-115
5.2.1 Cleaning and Checking the Air Filter.....	5-117
5.2.2 Removing the Air Filter.....	5-117
5.2.3 Cleaning the Air Filter.....	5-118
5.2.4 Daily Instrument Report	5-118
5.2.5 Cleaning / Decontamination	5-120
5.3 Module Replacement	5-121
5.3.1 General	5-121
5.3.2 Thermometer Removal	5-122
5.3.3 Reuniting Separated Liquid in the Thermometer	5-122

Table of Contents

6 – System Alert Messages

6.1 Error / Alert Messages	6–123
6.2 System Alerts (E error codes)	6–124
6.3 Workflow Alerts (W error codes)	6–137
6.4 Event Log Messages.....	6–150
6.4.1 LIS Related Messages	6–150

7 – Performance, Limitations, Organism Codes, and Panel Information

7.1 Performance Characteristics	7–161
7.1.1 Gram Negative Performance	7–161
7.1.2 Gram Positive Performance	7–164
7.1.3 Streptococci Performance (with BD Phoenix SMIC/ID, SMIC Panels)	7–166
7.1.4 Identification of Yeast Species	7–168
7.2 Limitations of the Procedure	7–168
7.3 Organism Codes, Short Names, Long Names.....	7–169
7.4 List of Antimicrobial Agents in BD Phoenix Panels	7–193
7.4.1 Gram Negative	7–193
7.5 List of Reagents and Principles Employed in the BD Phoenix System.....	7–200
7.5.1 Gram Negative	7–200
7.5.2 Gram Positive.....	7–202
7.5.3 Streptococci Panel	7–203
7.5.4 Yeast Panel.....	7–204
7.6 Taxa for ID/AST Determination.....	7–206
7.6.1 Gram Negative (0.5 McFarland).....	7–206
7.6.2 Gram Negative (0.25 McFarland).....	7–210
7.6.3 Gram Positive (0.5 McFarland)	7–212
7.6.4 Gram Positive (0.25 McFarland)	7–215
7.6.5 Streptococci	7–217
7.6.6 Yeast	7–218

8 – Glossaries

8.1 Glossary	8–221
8.2 Supplemental ID Test Abbreviations.....	8–223
8.3 Symbol Glossary	8–224

Table of Contents

9 – Replacement Parts

10 – Software Update Log

11 – International Contacts

1 – Introduction

1.1 Intended Use

The BD Phoenix™ Automated Microbiology System is intended for the rapid identification (ID) and antimicrobial susceptibility testing (AST) of clinically significant bacteria. The BD Phoenix system provides rapid results for most aerobic and facultative anaerobic gram-positive bacteria as well as most aerobic and facultative anaerobic gram-negative bacteria of human origin. The BD Phoenix System is also intended for the rapid identification of yeast and yeast-like organisms.

1.2 Summary and Explanation of the Test

Micromethods for the biochemical identification of microorganisms were reported as early as 1918.¹ Several publications reported on the use of the reagent-impregnated paper discs and micro-tube methods for differentiating enteric bacteria.¹⁻⁹ The interest in miniaturized identification systems led to the introduction of several commercial systems in the late 1960s, and they provided advantages in requiring little storage space, extended shelf life, standardized quality control, and ease of use. Many of the tests used in the BD Phoenix ID panels are modifications of the classical methods. These include tests for fermentation, oxidation, degradation and hydrolysis of various substrates. In addition to these, the BD Phoenix System utilizes chromogenic and fluorogenic substrates as well as single carbon source substrates in the identification of organisms.^{10, 11, 26-28}

The modern broth microdilution test used today has origins in the tube dilution test used in 1942 by Rammelkamp and Maxon to determine *in vitro* antimicrobial susceptibility testing of bacterial isolates from clinical specimens.¹² The broth dilution technique involves exposing bacteria to decreasing concentrations of antimicrobial agents in liquid media by serial two-fold dilution. The lowest concentration of an antimicrobial agent in which no visible growth occurs is defined as the minimal inhibitory concentration (MIC).

The introduction in 1956 of a microtitrator system, using calibrated precision spiral wire loops and droppers for making accurate dilutions, rapidly allowed Marymont and Wentz to develop a serial dilution AST test.¹³ The microtitrator system was accurate and allowed the reduction in volumes of antimicrobial agents. The term microdilution appeared in 1970 to describe the MIC tests performed in volumes of 0.1 mL or less of antimicrobial solution.¹⁴

The BD Phoenix AST is a modified miniaturized version of the micro-broth doubling dilution technique. Susceptibility testing in the BD Phoenix System is performed through determination of bacterial growth in the presence of various concentrations of the antimicrobial agent tested with the aid of the AST indicator in continuously incubated and read micro-wells in the BD Phoenix panels.

1.3 Principles of the Procedure

A maximum of 50 identification and antimicrobial susceptibility tests can be performed in the BD Phoenix M50 instrument at a time using BD Phoenix combination panels. A sealed and self-inoculating molded polystyrene tray with 136 micro-wells containing dried reagents, serves as the BD Phoenix disposable. The combination panel includes an ID side with dried substrates for bacterial or yeast identification and an AST side with varying concentrations of antimicrobial agents, growth and fluorescent controls at appropriate well locations. The BD Phoenix system utilizes an optimized colorimetric redox indicator for AST, and a variety of colorimetric and fluorometric indicators for ID. The AST broth is cation-adjusted (e.g., Ca⁺⁺ and Mg⁺⁺) to optimize susceptibility testing performance.

The BD Phoenix Panel is comprised of a 51 well ID side and an 85 well AST side. The ID side contains 45 wells with dried biochemical substrates and 2 fluorescent control wells. The AST side potentially contains up to 84 wells with dried antimicrobial agents and 1 growth control well. Panels are available as Emerge, ID only, AST only, or ID/AST combination. Unused wells are reserved for future use. BD Phoenix panels are inoculated with a standardized inoculum. Organism suspensions must be prepared only with the BD BBL™ CrystalSpec™ Nephelometer, the BD PhoenixSpec™ Nephelometer, or the BD Phoenix AP instrument. Once inoculated, panels are placed into the instrument and continuously incubated at 35 °C. The instrument tests panels every 20 minutes: on the hour, at 20 minutes past the hour, and again at 40 minutes past the hour for up to 16 hours if necessary. BD Phoenix panels are read only by the BD Phoenix M50 instrument. BD Phoenix panels cannot be read manually.

Organism Identification

The ID portion of the BD Phoenix panel utilizes a series of conventional, chromogenic, and fluorogenic biochemical tests to determine the identification of the organism. Both growth-based and enzymatic substrates are employed to cover the different types of reactivity within the range of taxa. The tests are based on microbial utilization and degradation of specific substrates detected by various indicator systems. Acid production is indicated by a change in phenol red indicator when an isolate is able to utilize a carbohydrate substrate. Chromogenic substrates produce a yellow color upon enzymatic hydrolysis of either *p*-Nitrophenyl or *p*-Nitroanilide compounds. Enzymatic hydrolysis of fluorogenic substrates results in the release of a fluorescent coumarin derivative. Organisms that utilize a specific carbon source reduce the resazurin-based indicator. In addition, there are other tests that detect the ability of an organism to hydrolyze, degrade, reduce, or otherwise utilize a substrate.

A complete list of taxa that comprises the BD Phoenix Database is provided in Section 7.6 Taxa for ID/AST Determination. Reactions employed by various substrates and a brief explanation of the principles employed in the BD Phoenix Gram Negative, Gram Positive, *Streptococcus*, and Yeast ID reactions are described in Section 7 – Performance, Limitations, Organism Codes, and Panel Information.

Antimicrobial Susceptibility Testing

The BD Phoenix AST method is a broth based microdilution test. The BD Phoenix system utilizes a redox indicator for the detection of organism growth in the presence of an antimicrobial agent.¹⁵ Continuous measurements of changes to the indicator as well as bacterial turbidity are used in the determination of bacterial growth. Each AST panel configuration contains several antimicrobial agents with a wide range of two-fold doubling dilution concentrations. Organism identification is used in the interpretation of the MIC values of each antimicrobial agent. A complete list of taxa for which the BD Phoenix panels can provide AST results is shown in Section 7 – Performance, Limitations, Organism Codes, and Panel Information. The list of antimicrobial agents and concentrations available for susceptibility testing in the BD Phoenix system is provided in Section 7.4 List of Antimicrobial Agents in BD Phoenix Panels.

The components required for testing include:

- 1 BD Phoenix **panels with panel closures**
- 2 BD Phoenix **ID Broth**
- 3 BD Phoenix **AST Broth**
- 4 BD Phoenix **AST Indicator Solution**
- 5 BD Phoenix AST-S Broth (for use with BD Phoenix Strep panels only)
- 6 BD Phoenix AST-S Indicator Solution (for use with BD Phoenix Strep panels only)
- 7 BD Phoenix **Inoculation Station**
- 8 BD Phoenix **Panel Caddy**

-
- 9 BD BBL CrystalSpec Nephelometer, the BD PhoenixSpec Nephelometer, or the BD Phoenix AP instrument
 - 10 25 µL pipettor and sterile tips
 - 11 Miscellaneous lab supplies (listed under Materials Required but Not Provided: in Section 4.3 Preparing Panels).

The BD Phoenix panel is for use with the BD Phoenix M50 instrument for either organism identification (ID) or antimicrobial susceptibility testing (AST), or for the combination of both. The panel is inclined with the inoculation ports at the top for filling. Separate inocula are added manually to the ID and AST ports. The inocula flow down the panel in serpentine fashion, filling the panel wells as the liquid front progresses toward the pad. The pad absorbs excess inoculum. Closures are manually inserted in the fill ports. An air admittance port is located in the divider area of the panel lid to ensure adequate oxygen tension in the panel for the duration of the test.

Principles of BD Phoenix AST Tests for the Detection of Resistance Markers

The following sections outline the principles of the BD Phoenix AST System in the detection of resistance markers in gram-negative or gram-positive organisms, including 1) detection of ESBL production among species of *Enterobacteriaceae*; 2) detection of vancomycin resistance in *Enterococcus* species (VRE); 3) detection of high-level aminoglycoside resistance in *Enterococcus* and *Streptococcus* species (HLAR); 4) detection of methicillin-resistance in staphylococci (MRS); 5) detection of β -lactamase production in *Staphylococcus* species (BL); 6) detection of macrolide resistance in *Streptococcus* species (MLSb); 7) detection of *mecA*-mediated Resistance with *S. aureus* (*mecA*); 8) detection of Vancomycin Resistant *Staphylococcus aureus* (VRSA); 9) detection of BD Phoenix Inducible Macrolide Resistance (IMLS) in *Staphylococcus* spp. For further information, consult the BDxpert manual.

BD Phoenix Extended Spectrum β -Lactamase (ESBL) Test¹⁶

The BD Phoenix ESBL test evolved from published data of known ESBL antibiogram patterns in the current literature.¹⁸⁻²¹ Selected strains of various species with known β -lactamase genotype/phenotypes in the family *Enterobacteriaceae*, including *Escherichia coli*, *Klebsiella* species (spp.), *Citrobacter* spp., *Enterobacter* spp., *Proteus* spp., and *Serratia* spp., were used to develop the BD Phoenix ESBL test. The BD Phoenix ESBL test is based on the principle of a differential response between the inhibitory effect of selected second or third generation cephalosporins in the presence or absence of a β -lactamase inhibitor, clavulanic acid. The principles of BD Phoenix ESBL test is similar to the CLSI ESBL broth microdilution confirmatory test.²² The BD Phoenix ESBL test is applied to *E. coli*, *K. pneumoniae* and *K. oxytoca*. Additionally, at the customers discretion, it can be applied to other enteric species where production of ESBL has been reported in literature. When a test result of ESBL is positive, the categorical interpretation of all penicillins, cephalosporins (except cephamycins), and aztreonam on the same BD Phoenix panel will be changed to R with BDxpert rule 1529. Customers can enable specific rules to be reported as tested.²²

BD Phoenix Vancomycin Resistant Enterococci (VRE) Test

The BD Phoenix VRE test is based on the SIR interpretation of vancomycin. The breakpoint selected in the instrument configuration is used for the categorical interpretation. The BD Phoenix VRE test was developed and optimized to match the CLSI standard broth microdilution test.^{22, 23} Selection of a breakpoint other than CLSI may result in less than optimal performance due to differences in categorical interpretations. Only *Enterococcus faecalis* and *E. faecium* with acquired resistance (*vanA* or *vanB*) will be reported as positive.²²

BD Phoenix High-Level Aminoglycoside Resistance (HLAR) Tests

The BD Phoenix HLAR tests for *Enterococcus* are based on the growth response in a single well containing either a high-level concentration of gentamicin or streptomycin. These tests were developed and optimized against both the CLSI standard broth microdilution and the CLSI screening agar test.²²

The BD Phoenix HLAR tests for *Streptococcus* are based on the growth response in a single well containing gentamicin, kanamycin, or streptomycin. These tests were developed and optimized using the CLSI recommended standard broth microdilution.

BD Phoenix Methicillin-Resistance in Staphylococci (MRS) Test

The BD Phoenix MRS test is based on the SIR interpretation of oxacillin with *Staphylococcus* species. When an MRS test result is positive, several BDxpert rules are designed to handle the reporting and the interpretations of all beta-lactam drugs. BD Phoenix cefoxitin MIC result is used to predict *mecA*-mediated resistance in *Staphylococcus aureus*. A special BDxpert rule is designed to report MRS using cefoxitin results for *Staphylococcus aureus*. The surrogate drug, cefoxitin, has been validated as a better indicator for the presence of *mecA* in staphylococci. The breakpoint selected in the instrument configuration is used for the categorical interpretation.

BD Phoenix Gram-Positive β -lactamase (BL) Test¹⁶

The BL test available in the BD Phoenix AST System is a nitrocefin based β -lactamase test. The nitrocefin based test is a direct detection method located on the ID side of the BD Phoenix panel. The performance of this test was established against the results of testing with BD BBL Cefinase™ Discs (Cat. No. 231650) as the reference method. Currently, only *Staphylococcus* species will be evaluated with these tests. When the result of BL test is positive, the categorical interpretation of all penicillinase labile penicillins on the same Phoenix panels will be changed to resistant.²²

BD Phoenix Macrolide Resistance in Streptococci (MLSb) Test

The BD Phoenix Macrolide Resistance test is based on SIR interpretation of erythromycin and clindamycin. The breakpoint selected in the instrument configuration is used for the categorical interpretation. Erythromycin resistant and clindamycin resistant *Streptococcus* isolates will be reported as macrolide/lincosamide/streptogramin B (MLSb) phenotype.

BD Phoenix *mecA*-mediated Resistance Marker for *Staphylococcus aureus* (*mecA*)

The BD Phoenix *mecA* test is used to predict *mecA*-mediated resistance in *Staphylococcus aureus*. The principle is similar to the CLSI-recommended Disk Diffusion test, which uses a cefoxitin (FOX) disk to predict *mecA*-mediated resistance in *S. aureus*. The performance of the test was established against multiplex PCR methods²⁵ as well as the Disk Diffusion test. With the BD Phoenix *mecA* test, the *mecA*-specific FOX MICs used for detection of the resistance marker will be configured in the instrument. When the *mecA* resistance marker is detected, the interpretations for all beta-lactam drugs on the same BD Phoenix panel are changed to resistant,²² and the BD Phoenix *mecA* resistance marker is set.

BD Phoenix Vancomycin Resistant *Staphylococcus aureus* (VRSA) Test

The BD Phoenix VRSA detection is based on the SIR interpretation of vancomycin when testing *Staphylococcus aureus*. The breakpoint selected in the instrument configuration is used for the categorical interpretation. The BD Phoenix VRSA test was developed and optimized to match the CLSI standard broth microdilution test, and verified with known VRSA isolates. Selection of a breakpoint other than those found in CLSI M100-S25 may result in less than optimal performance due to differences in categorical interpretations. Only *Staphylococcus aureus* with true resistance (isolates

containing resistance marker such as *vanA* gene) will be reported as VRSA. Strains of *S. aureus* with vancomycin intermediate results (GISA/VISA) will be identified and reported by separate BDxpert rules. The BD Phoenix Gram Positive AST panel detected vancomycin resistance in the VRSA *S. aureus* strains available at the time of comparative testing. The ability to detect resistance in other *S. aureus* strains is unknown due to the limited number of resistant strains available for comparative testing.

BD Phoenix Inducible Macrolide Resistance (IMLS) Test in *Staphylococcus* species

The BD Phoenix Inducible Macrolide Resistance (IMLS) Test is used to detect inducible macrolide lincosamide-streptogramin B (MLSb) resistance in *Staphylococcus* species. MLSb resistance, usually encoded by *ermA* or *ermC* genes, may be either constitutive (always expressed) or inducible after exposure to a macrolide antibiotic (e.g. erythromycin, clarithromycin, etc.). The BD Phoenix Inducible Macrolide Resistance Test is based on the same principle as the CLSI recommended Disk Approximation Test (D-Test) for the detection of inducible clindamycin resistance. When the BD Phoenix Inducible Macrolide Resistance Test result is positive, the categorical interpretation of clindamycin on the same BD Phoenix panel will be reported as resistant and accompanied by a separate BDxpert message. *Staphylococcus* isolates resistant to both erythromycin and clindamycin on initial testing will be reported as constitutive MLSb resistance to distinguish them from isolates that are resistant to macrolides alone by efflux mechanism.

1.4 Precautions

For *in vitro* Diagnostic Use.

All patient specimens and microbial cultures are potentially infectious and should be treated with universal precautions. Please refer to CDC manual *Biosafety in Microbiological and Biomedical Laboratories*, 5th Edition, 2009, as well as other recommended literature.

Panels once inoculated should be handled carefully until placed in the instrument.

1.5 System Overview

1.5.1 Instrument Overview

1.5.1.1 Instrument Optical System and Drive Mechanism

Within the main instrument incubation bay, the cylindrical carousel and its drive are the only moving parts. The carousel is normally driven at either 1.0 RPM or 2.0 RPM. The drive speed, acceleration and deceleration are controlled by the central processor. The carousel is divided vertically into two tiers, each of which functions as a relatively independent optical source and detection system for panels placed on that tier. Each tier has its own microcontroller to control data acquisition and transmission. The two microcontrollers communicate with the central processor over a serial communications link.

Visible illumination in the red, green and blue spectral regions for each tier is provided by a Light Emitting Diode (LED) source board. The LED source currents are programmable in order to compensate for signal loss at the panel extremes due to parallax and other factors. A source monitor system averages signal from two onboard photodiodes to monitor visible source output. Fluorescent UV (ultraviolet) LEDs provide UV illumination.

1.5.1.2 Carousel Assembly

The Carousel Assembly is a cage-like structure comprised of aluminum rings and vertical ribs bolted together to form a right circular cylinder. The carousel holds 52 panel holders in two tiers. Each tier holds one normalizer panel and 25 sample panels.

During normal operation, the carousel is driven at either 1.0 RPM or 2.0 RPM, depending on the current operation. During panel location, the carousel can rotate at up to 10.0 RPM. The carousel rotates counterclockwise. A single rotation of the carousel is used to perform an inventory scan to identify panels within the instrument following a door closure. A complete test cycle requires 7 minutes.

1.5.1.3 Incubation System

The Incubation System maintains the carousel, panel carriers, and panels at a constant nominal temperature of 35 °C. The system is a recirculating forced-air convection design. There is a single filter to remove dust from the electronics bay. The system consists of a squirrel cage blower powered by a brushless DC motor, a coil-style heater section with automatic over-temperature shut-off, thermoformed inlet and return ductwork, and a user-replaceable polyester fiber air filter.

1.5.1.4 Panel Status and Internal Barcode Scanner Assembly

Panel status indication and panel identification is accomplished on the same tower in the instrument, mounted directly behind the carousel in the door area. Panel Status is indicated by red, green, and amber LEDs mounted on a panel behind the carousel that shine through the light pipes and panels. Two columns of the carousel are exposed when the instrument door is open.

Two barcode scanners capable of reading panel barcodes are located on the same tower as the indicators. Each scanner reads panel barcodes on its respective tier of the carousel. Panel barcodes are affixed at the top of the panel on the base (bottom) side of the panel, facing inward on the carousel when placed in the instrument.

1.5.1.5 External Barcode Scanner

The external barcode scanner can be used to read barcoded accession numbers that have been placed on the panels, as well as the panel's own sequence number barcodes. The accession barcodes can be used to link specimen identification information to specific panels in the instrument.

1.5.2 Software and Operation Overview


Figure 1-1 – BD Phoenix M50 Instrument

The All-In-One (AIO) PC presents all the information needed to monitor instrument status, to enter and remove panels, set up the instrument, print reports, and perform routine instrument maintenance. The AIO PC controls up to two BD Phoenix M50 instruments (see Figure 1-3 – BD Phoenix M50 User Interface for Two Instruments). The AIO PC provides connectivity to the hand-held barcode scanner and printer. The Interface to a BD Data Management system or LIS is optionally supported by the AIO PC.


Figure 1-2 – BD Phoenix M50 User Interface for One Instrument

The User Interface on the AIO PC is divided into 3 regions (see Figure 1-2 – BD Phoenix M50 User Interface for One Instrument). The application header and footer appear on all screens of the BD Phoenix M50 instrument user interface. The middle region corresponds to the selected tab in the footer; in the figure above the status tab has been selected. The Application Header contains status information for the BD Phoenix M50 instrument. Either one or two instrument summary groups are displayed in the center of the header (See Figure 1-3 – BD Phoenix M50 User Interface for Two Instruments). The application footer contains the tabs to navigate between different screens.


Figure 1-3 – BD Phoenix M50 User Interface for Two Instruments

1.5.3 Panel Overview

The BD Phoenix panel is available in three formats: ID only, AST only, and combination ID/AST testing. The pour and seal serpentine design is optimized for safety and leak-resistance. Each well in the disposable contains approximately 50 μ L of inoculum in an environment that prevents significant evaporation during the course of incubation.


Figure 1-4 – Panel Inoculation Station

A Panel inoculation station holds six tubes of broth (ID, AST) and three panels held at an angle of 24° in order to provide proper gravity-driven inoculum flow through the panel. See Figure 1-4 – Panel Inoculation Station.

A Panel Transportation Caddy is a molded plastic tray used to transport filled and sealed Phoenix panels from the preparation bench to the BD Phoenix M50 instrument. The caddy capacity is 20 panels. See Figure 1-5 – Panel Transportation Caddy.

Figure 1-5 – Panel Transportation Caddy

1.5.4 Testing Overview

After a door-opening/-closing event, the instrument reads panel barcode labels and then performs a scan using the red LEDs to determine if panels are present, and to locate or map well positions. Panel readings are made on the hour and at 20 and 40 minutes after the hour. The instrument test sequence begins with the system checking that the door is locked and the current time is read. Tiers perform dark readings, then the UV LEDs are turned on and allowed to warm up. UV readings are then taken for one revolution. The UV LEDs are turned off. Next the red LEDs are turned on, allowed to warm up, and red LED readings are taken. Then the green LEDs go through the same sequence. Finally, the blue LEDs go through this sequence.

If the test cycle completes successfully, the time is saved. A successful test cycle occurs when there are no carousel errors and the user does not preempt the test by initiating a Load Panels, Unload Panels, or Locate Panel operation. When the user initiates these operations, the current time is compared to the last test cycle's start time. If more than 30 minutes have not elapsed, the requested operation (for example: Load Panels, Unload Panels) is performed. If more than 30 minutes have elapsed since the last test cycle's start time, a test cycle must successfully complete before the user is permitted to perform a Load Panels, Unload Panels, or Locate Panel operation. The Panel In/Out indicator is off (see Figure 2-7 – Panel In/Out Indicator) and the **cannot get into instrument** tone sounds when a panel operation is requested.

After each test:

- The summary counters on the Status screen are adjusted to indicate current statuses
- Panels/records that require user action have a Needs Attention set
- System alerts are reported in the System Alerts list
- Auto Association occurs

The screen displays the icon below when a panel test is taking place. The number next to the icon indicates the number of minutes remaining in the test cycle.

1.5.5 Normalizers

Normalizers serve as reference panels for adjustment of the instrument's optical detection system and they are used for adjustment of the LEDs (red, green, blue and UV). For red, green, blue, and UV correction, data from the normalizer panels is used to correct for variations in optical channel gain and to compensate for well-to-well parallax. Raw signals are source monitor corrected and are then ratioed to the value of the corresponding well of the normalizer panel. This ratio is multiplied by a correction factor, which is the expected transmittance level of the normalizer, to scale the resultant values. The fluorescence (UV) signal from the normalizer is used to ensure that the UV intensity reaching the panel is within the acceptable range for proper interrogation of the fluorescent ID substrates on the BD Phoenix panel.

Normalizers expire two years after installation. A system alert will occur before expiration to allow for scheduling replacement.

The normalizer panels are always located in station number 0 of each tier. When the door is opened, routine carousel access does not present the normalizer panel. The normalizer panel is constrained in its carrier to prevent inadvertent removal.

Automatic Adjustment of Light Sources

Automatic light source adjustment attempts to bring normalizer readings within acceptable ranges. It is performed when the following conditions are met: the instrument is warmed up, idle, and there are no ongoing panels. UV adjustment prohibits access to the instrument for activities such as entering panels into unaffected tiers, removing panels, and performing maintenance checks.

Based on the results of the system's Built-In-Test (BIT), as well as time and power cycling factors, the system detects that readings from the ultraviolet (UV) and visible light sources (testing LEDs) are out of tolerance. There are two main levels to this condition:

- 1 The deviation is great enough that panel results are invalidated (panel testing for the tier is aborted); the tier's stations are blocked; an automatic light source adjustment is performed as soon as all panels in the instrument complete testing.
- 2 The deviation is within limits that do not affect panel results; however available stations are blocked; as soon as all panels on the tier complete testing or are removed, an automatic light source adjustment is performed.

If any of the above conditions occur, E type error codes, and sub-codes will provide information for each of the tiers.

1.6 Limited Warranty

This warranty details specific legal rights. Additionally, there may be other rights that vary from region to region.

The BD Phoenix M50 Automated Microbiology System is warranted to the original purchaser to be free from defects in materials and workmanship for a period of one year following installation. BD's sole responsibility under this warranty shall be to repair or replace any instrument or its components (except for expendable supplies such as printer cartridges, paper, or filters) which under normal operating conditions, prove to be defective within one year of delivery.

BD will furnish new or re-manufactured components upon its option. All replacements shall meet new part specifications and shall be warranted as above for the remainder of the one year period. Replaced components become the property of BD.

It is understood that the equipment covered by this Agreement has been installed in accordance with the recommendations and instructions in the BD Phoenix M50 System User's Manual.

Any damage to a BD Phoenix M50 instrument resulting from the insertion or removal of cables that connect this instrument to systems other than those approved or supplied by BD, or the failure of the owner to maintain reasonable care and precautions in the operation and maintenance of the system, will void this warranty and terminate the obligations of the manufacturer as stated herein.

This warranty is in lieu of all other warranties, whether express or implied, including but not limited to, warranties of merchantability, or fitness for a particular use. In no event will BD be liable for indirect, incidental, special or consequential damages regardless of whether BD has been advised of such.

1.7 Conventions

1.7.1 Notes, Cautions, and Warnings


Throughout this manual, important information is presented in boxes offset from the regular text, and is labeled as either a NOTE, CAUTION, or WARNING. These messages are formatted as shown below and bear the following significance:

NOTE

Important information about system use worthy of special attention is presented as a NOTE.

CAUTION

Information on an activity which potentially could cause damage to the instrument or system is presented as a CAUTION.

1.7.2 Summary of Cautions and Warnings

- Protection provided by this equipment may be impaired if the equipment is used in a manner not consistent with the instructions in this manual.
- Due to the size and weight of the BD Phoenix M50 instrument, two people are required to lift the instrument in the absence of mechanical lifting devices.
- All system users should become thoroughly familiar with all controls and indicators before attempting to operate the instrument.
- Observe established precautions against microbiological hazards throughout all procedures. All specimens should be handled according to CDC-NIH recommendations, CLSI guidelines, or local institution guidelines for any potentially infectious human serum, blood, or other body fluids. Prior to discarding, sterilize specimen containers and other contaminated materials by autoclaving.
- In addition to wearing gloves, the use of disposable lab coats or gowns and protective glasses or goggles is recommended when working around the instrument.
- The instrument door is electromagnetically latched and is controlled by the instrument software. Never attempt to defeat the door latching mechanism, or to open the door when the unlocked icon is not displayed. Serious injury can be caused by the rotating carousel manually. If the carousel is not completely stopped when the door is opened, immediately contact BD for service. Never attempt to rotate the carousel manually or serious injury may result.
- When the system displays alerts and errors, immediately respond to the condition.
- All maintenance and repair other than the procedures described in Section 5.2 Routine Maintenance must be performed by qualified service personnel. Non-compliance with this warning may result in personal injury or instrument malfunction.
- All portions of the body that could possibly come into contact with the affected instrument surfaces must be completely covered before beginning the decontamination process.
- If any error sub-codes other than those listed here appear, note the sub-code and contact BD for assistance.
- If the recommended corrective actions do not solve the problem, contact BD.

1.8 Bibliography

- 1 Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. *Am. J. Public. Health.* 8:922-923.
- 2 Arnold, W.M., Jr., and R.H. Weaver. 1948. Quick microtechniques for the identification of cultures. I. Indole production. *J. Lab. Clin. Med.* 33:195-195.
- 3 Bachmann, B., and R.H. Weaver. 1951. Rapid microtechniques for identification of cultures. V. Reduction of nitrates. *Am. J. Clin. Pathol.* 21:195-196.

-
- 4 Hannan, J., and R.H. Weaver. 1949. Quick microtechniques for the identification of cultures II. Germentations. J. Lab. Clin. Med. 33:1338-1341.
 - 5 Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
 - 6 Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
 - 7 Snyder, M.L. 1954. Paper discs containing entire culture medium for the differentiation of bacteria. Pathol. Bacteriol. 67:217-226.
 - 8 Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Publ. Hlth. Trop. Med. 25:96-100.
 - 9 Weaver, R.H. 1954. Quicker bacteriological results. Am. J. Med. Technol. 20:14-26.
 - 10 Kampfer, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877-2879.
 - 11 Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
 - 12 Rammelkamp, C.H., and T. Maxon. 1942. Resistance of *Staphylococcus aureus* to the action of penicillin. Proc. Soc. Biol. And Med. 51:386-389.
 - 13 Marymont, J.H. and R.M. Wentz. 1966. Serial dilution antibiotic sensitivity testing with the microtiterator system. Am. J. Clin. Pathol. 45:548-551.
 - 14 Gavan, R.L., and M.A. Town. 1970. A microdilution method for antibiotic susceptibility testing: an evaluation. Am. J. Clin. Pathol. 53:880-885.
 - 15 Lancaster, M.V. and D. Rebecca. 1996. Antibiotic and cytotoxic drug susceptibility assays using resazurin and poisoning agents. US Patent #5,501,959.
 - 16 Data on file at Becton, Dickinson and Company.
 - 17 Jorgensen, James H., et al. ed., Manual of Clinical Microbiology, 11th Edition, ASM Press, Washington, D.C., 2015.
 - 18 DM Livermore. 1995. Beta-Lactamases in laboratory and clinical resistance. Clin. Microbiol. Rev. 8: 557-584.
 - 19 K Bush, GA Jacoby, and AA Medeiros. 1995. A functional classification scheme for beta-lactamases and its correlation with molecular structure. Antimicrob. Agents Chemother. 39: 1211-1233.
 - 20 GA Jacoby and P Han. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. J. Clin. Microbiol. 34: 908-911.
 - 21 Kenneth S. Thomson, Christine C. Sanders, and Ellen Smith Moland. 1999. Use of Microdilution Panels with and without Beta-Lactamase Inhibitors as a Phenotypic Test for Beta- Lactamase Production among *Escherichia coli*, *Klebsiella* spp., *Enterobacter* spp., *Citrobacter freundii*, and *Serratia marcescens*. Antimicrob. Agents Chemother. 43: 1393-1400.
 - 22 CLSI. M100-S26 Performance Standards for Antimicrobial Susceptibility Testing; Twenty-sixth Informational Supplement. January 2016.
 - 23 CLSI. M07-A10 Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard – Tenth Edition. January 2015.
 - 24 Livermore, D.M., T.G. Winstanley, and K.P. Shannon. 2001. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. J. Antimicrob. Chemother. 48, Suppl. S1, 87-102.
 - 25 Cauwelier, et al. "Evaluation of a disk diffusion method with cefoxitin (30 mcg) for detection of methicillin-resistant *Staphylococcus aureus*." Eur J. Clin. Microbiol Infect Dis. May 2004,23(5):389-392.
 - 26 Barnett, J.A., Payne, R.W. and Yarrow, D., Yeasts: Characteristics and Identification, 3rd Edition, Cambridge University Press, Cambridge, U.K. 2000.
 - 27 Larone, D. H., Medically Important Fungi, A Guide to Identification, 5th Edition, ASM Press, Washington, D.C. 2011.
 - 28 de Hoog, G.S., Guarro, J., Gené, J. and Figueras, M.J., Atlas of Clinical Fungi, 4th Edition, Centraalbureau voor Schimmelcultures, Utrecht, the Netherlands/Universita Rovira I Virgili, Reus, Spain, 2014.
-

2 – Controls, Icons, Indicators, and Installation

2.1 General

This section describes instrument specification, installation, software icons, controls and indicators, and user configuration of the BD Phoenix M50 instrument.

WARNING

PROTECTION PROVIDED BY THIS EQUIPMENT MAY BE IMPAIRED IF THE EQUIPMENT IS USED IN A MANNER NOT CONSISTENT WITH THE INSTRUCTIONS IN THIS MANUAL.

CAUTION

The intake filter at the lower front right corner of the BD Phoenix M50 instrument must remain unobstructed at all times. Restricted air flow may cause excessive temperatures in the instrument, which can affect test results and possibly cause hardware malfunctions.

2.2 System Startup

Whenever power is applied to the system, it is initialized, performs self-diagnostics, and reports any problems to the user through alerts. If any files are missing or corrupted which would prevent proper operation of the system, the startup process is aborted. If not, the computer loads the operating system and user interface, and begins the warmup period (indicated by the Instrument is warming up icon). Afterwards, the system awaits the initiation of panel testing.

The temperature standard panel should be left in the instrument at least 45 minutes before reading it when the instrument is first powered up.

2.2.1 Instrument Installation

The BD Phoenix M50 Instrument is to be installed only by BD representatives.

CAUTION

Due to the size and weight of the BD Phoenix M50 instrument, at least two people are required to lift the instrument in the absence of a mechanical lifting device.

Language

The language setting is performed at installation by a BD representative. Supported languages include:

- English
- Spanish
- French
- Portuguese
- German
- Japanese
- Italian

The BD Phoenix M50 instrument should be installed in an area that is free from undue vibration, direct sunlight, high humidity, dust, temperature extremes, and corrosive or explosive vapors or gases. The system will operate within specifications in room temperatures between 18–30 °C (64.4–86 °F). Relative humidity should be between 20% and 80% (non-condensing). The left, rear, and right sides of the instrument should be placed at least three inches from any wall. Environments which exceed these limits could adversely effect the performance of the system components.

The carousel should maintain its temperature to within plus or minus 1.5 °C of the temperature controller's setting (35 °C). This accuracy can be assured only if the room temperature meets the requirements given above.

Use of earthquake anchoring is strongly recommended in areas susceptible to earthquake activity.

A stacking kit is required when stacking two instruments on top of one another.

Installation Category II and Pollution Degree 2 as per IEC 664.

2.3 Instrument Specifications

2.3.1 Specification Chart

Physical Dimensions	
Height	21 in (535 mm)
Width	32 in (815 mm)
Depth	30 in (765 mm)
Clearance–Right	21.5 in (546 mm)
Clearance–Left, Rear	3 in (76mm)
Clearance–Front	18 in (457 mm)
Weight–Empty	118 lb (53.5 kg)
Weight–Full	127 lb (57.6 kg)

Environmental Requirements	
Non-Operating Storage	
Temperature	-17.8–65 °C (0–149 °F)
Humidity	20–80% RH, non-condensing
Operating Conditions	
Temperature	18–30 °C (64.4–86 °F)
Humidity	20–80% RH, non-condensing
Locations	Level surface, No direct sunlight, No direct heat
Altitude	2000 meters
Noise	at 1 m ≤ 55 dBA using ANSI Type 2 sound meter

Electrical Requirements	
Input Voltage	100–240 VAC
Input Current	6 Amp
Input Line Frequency	50/60 Hz

Instrument Heat Output at Input Voltage = 120 VAC	
Ambient Temperature	BTU/hr
18 °C	648
24 °C	464
30 °C	287

Optical Specifications	
Peak Wavelengths	
428 nm to 623 nm, visible spectrum	
375 nm, UV Excitation	
410 nm -- 640 nm Bandpass, UV Emission	

2.3.2 Controls and Indicators

2.3.2.1 Power Switch and Inputs / Outputs

Disconnecting power from the instrument:

- 1 Place power switch in the **OFF** position by pressing switch on the O side.
- 2 Unplug power cord from the instrument.


Figure 2-1 – Power off and disconnected

2.3.2.2 Connector Plate


Figure 2-2 – BD Phoenix M50 Connector Plate

USB Connector	Serial (LIS) Connector
USB Connector	Remote Alarm Connector
USB Host Connector (to All-In-One Computer)	Network Connector

2.3.2.3 All-in-one (AIO) PC controls


Figure 2-3 – All-in-one Controls

The AIO icons are located on the lower front portion of the PC. The control buttons can be found directly under each icon on the bottom portion of the standing PC.

2.3.2.4 Power Controls


Figure 2-4 – Power On and Power Off

2.3.2.5 Other Indicators

Figure 2-5 – Symbols on Body of BD Phoenix M50 Instrument

From left to right:

- Biohazard
- Electrical Recycling required
- On/Off

2.3.2.6 Instrument Alert Indicator


Figure 2-6 – Instrument Alert indicator

The instrument alert indicator represents the current alert status. See the information below for status details:

Instrument Alert Indicator		
LED Appearance	State	Meaning
Off	No Alerts / Instrument Communicating	The instrument does not have any outstanding alerts at this time and is connected to the AIO PC.
Blinking	Instrument not connected	The instrument is not connected or it is not communicating with the AIO PC.
Solid On	Instrument Alert	The instrument has outstanding alerts. Details are available via the alert screen. The instrument is connected to the AIO PC.

2.3.2.7 System Alerts

As activities are performed on the BD Phoenix M50 instrument, and as testing progresses, system alerts and errors may occur. Different types of alerts and errors are flagged by one or more of the following: **E** and **W** error codes (see Sections 6.2 System Alerts (E error codes) and 6.3 Workflow Alerts (W error codes), audible tones (see Section 2.3.2.8 Audible Tones and Alarms), the system Alert icon appearing on the LCD screen, or the instrument's System Alert indicator flashing. A System Alert icon (see Section 2.5 Software Icons) in the upper left corner of the screen indicates the presence of a system alert. Touch the icon to display the System Alerts screen. The System Alerts screen enables the review any existing system alerts that may have occurred or that may still exist in the instrument.

Touching the >> icon expands the alert message providing more detailed information. All the **E** type errors are listed in Section 6.1 Error / Alert Messages. Not all error sub-codes can be addressed by the user.

CAUTION

Immediately respond when the system displays alert and error notifications.

2.3.2.8 Audible Tones and Alarms

As discussed above, audible tones and alarms sound to inform the user of various operational states of the instrument. Refer to the chart below for a detailed description of such tones and alarms.

Type	Example	Sound
Informational		
Acknowledge	Scanning a barcode	Single short high beep
Sample audible alert	In Configuration mode, the audible alert volume was adjusted	Three tones progressing from high to low ("Figaro")
Carousel halted	The carousel has stopped, the door can be opened	Three tones progressing from low to high
Door closed	The door has been completely closed and latched	Two short high beeps
Activity completed	An operation has been completed	Three short high beeps

Type	Example	Sound
Alarm		
Activity error (volume configurable)	Incorrect barcode scanned	Short high beep then short low beep – sequence repeated four times
System alert (volume configurable)	Optical failure	Single medium beep – one second on, three seconds off, repeating
Error tone	The attempted action cannot be performed	Short high beep then short low beep
Door open alarm	Door has remained open longer than 5 minutes	Continuous shrill trill
Cannot get into instrument	A higher priority activity prevents user access to the instrument	Single low beep
Critical Panel/ Resistance Marker notification	Critical Panel: partial results, panel completes, or ID is determined; Resistance Marker: marker has been triggered	Tone 1 (Medium frequency): Long beep, short beep, long beep, short beep; Tones 2 (low frequency) and 3 (very high frequency): Long beep, short beep, long beep, long beep
Carousel jammed	Carousel motion is impeded	High pitched trill

2.3.2.9 Instrument Door/Panel Button Status


Figure 2-7 – Panel In/Out Indicator

The instrument has an indicator that provides the current panel button and door status to the user. The chart below explains the status details (refer to Section 3.3.2 Inserting Panels in the Instrument).

Panel In/Out Indicator		
LED Appearance	State	Meaning
Off	Panels In/Panels Out not available	If the instrument is powered: Panels In/ Panels Out are not available at this time. The door is locked and may not be opened by the user. The indicator is off when the instrument is performing an operation that may not be interrupted by a Panels In or Panels Out operation.
Blinking	The door on the instrument is unlocked	The door is unlocked and the user may open the door to access the panels.
Solid On	Panels In/Panels Out active	The Door is locked and may not be opened by the user. The user may press either the Panels In or Panels Out button.

2.3.2.10 Station Indicators

Each accessible station has a set of LED indicators providing the station or panel's status. The station indicators are located in the center of the station (see Figure 2-8 – Station Indicators). The appearance of green, amber, or red indicate the conditions shown in the table below for a given station.


Figure 2-8 – Station Indicators

Station Status Indicators	
LED Appearance	Meaning
Off	Ongoing Panel (if panel is in station) or Available Station (if station is empty)
Green	Removable Panel
Amber	Panel found using locate panel function
Red	Blocked Station (also Temperature Standard Panel)

2.3.3 Printer

For an explanation of controls and indicators on the printer, refer to the manufacturer's operating instructions which is furnished separately. Note that a local printer should be connected to one of the instrument's USB ports. One printer can serve two BD Phoenix M50 instruments connected via one AIO PC.

2.3.4 Barcode Scanners

A hand-operated scanner is provided with the system to facilitate scanning panel barcodes. To scan a barcode press the trigger on the bottom side of the scanner. A single beep indicates a successful scan. Press the trigger to scan each additional barcode.

2.3.5 Touchscreen Keyboard

A touchscreen keyboard is available on the AIO PC for typing data directly into screen fields. To access the touchscreen keyboard touch the field where data is to be entered. Then touch the keyboard icon that appears to expose a full size touchscreen keyboard. Data can be entered directly on the touchscreen or a regular keyboard can be used. When all data is entered via the touchscreen, touch the return key to close the touchscreen keyboard.

2.3.5.1 Enabling/Disabling Touchscreen Fields

The touchscreen user interface requires certain fields to be selected in order for them to be enabled or disabled. To enable or disable a field, simply touch a box and/or circle to make your selection. When the field is represented by a box, touching this box will insert a check indicating that the feature has been enabled. Touching this field again will remove the check; therefore indicating that the feature is disabled.

2.4 User Access and Management Summary

User Access

User login is accessed by touching **Log In** on the Status screen of the BD Phoenix M50 instrument application. When a user is logged in, the text changes to Log Out. The current user name is displayed in the upper right corner of application header region.

Logging In

Select **Log In** to access the user log in screen. The user name is the full user name. The tech ID will not be accepted as the user name for logging in. If improper credentials are submitted, the system will

initiate an alert and prompt for retry. After 5 unsuccessful login attempts, the user account will be locked. If the user landed at this screen inadvertently or wishes to exit without logging in, select **Cancel** from the main screen.

NOTE

If the user has logged in, but the screen is idle for 15 minutes, the user will be logged out and any unsaved data will be lost.

Logging out

To access the log out screen, select **Log Out** from the Status screen and a confirmation will appear verifying that the user wishes to log out. To remain logged in, select **Cancel** to navigate back to the main application screen.

User Role Summaries

- **Administrator** — Lab administrator for lab location. All system functionality is accessible including system configuration. Administrators can create and delete other lab administrators and general lab users.
- **General User** — All system functionality is accessible except system configuration, which is read only. Users can reset their own password.

User Management

Access the user management features by navigating to the Configuration tab. The Users tab is first in the Configuration screen tab list.

General Lab Users

The only management feature available to General Lab users is the ability to change their passwords. All users will be required to change their passwords at certain intervals. To change passwords enter the new password and enter it again in the **confirm** field. Then select **Update Password**.

Lab Administrator Users

Users that are at the Lab Administrator level will be able to access the full set of user management features.

A list of both active and inactive users of the system is displayed for Lab Administrators. User Details show the User ID Role, and Tech ID for the selected user.

Administrators can select a user from the list and perform one of the following actions:

- Reset a user's password.
- Deactivate/disable a user account. Deactivated User accounts still appear in the list.
- Add new users to the system.

2.5 Software Icons

* Icon does not appear when connected to BD EpiCenter™.

** Icon does not appear when in standalone mode

GENERAL				
Collapsed	Expanded	Enabled	Print	Help

STATUS SCREEN			
Log in/Log Out	Silence Alarm	Status	Panel Login
Results	*Finalization	Inventory	Reports

ALERTS/NOTIFICATIONS			STATUS STATION	
*Cannot Finalize	System Alert	Special Message	Blocked Station	Empty Station

PANELS				NEEDS ATTENTION			
Ongoing Panel	Removable Panels, Complete	Final ID Edited, AST Result Edited	*Resistance Marker	Rapid Complete	Active Alert	No Active Needs Attention	Active Needs Attention

INSTRUMENT STATUS					
 Door Locked	 Door Unlocked	 Idle	 Processing Data	 Warming Up	 Temperature
 Testing	 Instrument Number	 Light Source Adjustment (RBV, IUV)	 Inventory Scan	 Isolation Recovery (Uploading Data)	 **EpiCenter Configured and Connected
 **EpiCenter Configured but not Connected	 *LIS Configured and Connected	 *LIS Configured but Disconnected	 Instrument Connection Status – Connected (All-in-One PC to Instrument)	 Instrument Connection Status – Disconnected (All-in-One PC to Instrument)	

2.6 Software Setup

The system ships with all setup parameters preset to factory default values. However, before using the instrument for panel testing, review the setup parameters to see if they are suitable for the laboratory. These parameters are grouped in the following categories:

- User Management (Section 2.4 User Access and Management Summary)
- Instrument Configuration (Section 3.10.10 Instrument Sub-Tab)
- Communications Configuration (Section 3.10.3 Communications Sub-Tab)
- Custom Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab)
- BDxpert Rule Configuration (Section 3.10.7 BDxpert Rules Sub-Tab)
- Rapid Reporting Configuration (Section 3.10.8 Rapid Reporting Sub-Tab)
- Panel Lot Definition (Section 3.10.9 Panel Lot Definition Sub-Tab)

Instrument setup parameters are explained below. The Status screen, which is active immediately after instrument startup, is shown in Figure 1-2 – BD Phoenix M50 User Interface for One Instrument. The Status screen is the default when no activity has been initiated. This screen's tabs are covered in detail in Section 3 – User Interface Operations.

Software updates are user installable. Insert the BD provided USB key in the USB connector on the AIO PC. On the Maintenance tab select Task Category: Software, Task: Upgrade Execute. The software update will be recorded on the Daily Instrument Report for reference. See Section 3.9.1 Maintenance Hardware Functions for step by step instructions.

NOTE

Custom Breakpoints

If any of the standard interpretation rule sets have changed when the software is updated or when a new BD Phoenix Update Data (PUD) is installed, then new rule sets are installed into the system database. Since custom rule sets are based on standard rule sets, a custom interpretation rule set will be merged with new standard rules. If custom breakpoints are used, be sure to print out the Custom Breakpoint Difference Report (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab) after each system update or PUD install.

2.7 Isolation Mode

Isolation Mode is the condition that exists when communication between the BD Phoenix M50 instrument and the AIO PC is lost. Isolation Mode is designed to allow the ID/AST system to avoid test cycle gaps when the AIO PC is temporarily disconnected from the instrument. During Isolation Mode, the instrument continues to collect panel test cycle data and incubate the panels. However, Isolation Mode is not intended to enable routine workflows such as Panel Login or Removing Completed Panels. Since panel result analysis occurs at the AIO PC, no panel results are available for completed panels while the system is in Isolation Mode.

Please note the following conditions about Isolation Mode related to system operation:

- The yellow instrument alert indicator will blink when the instrument is in Isolation Mode.
- The AIO PC displays errors when communication is lost with the instrument. On the status screen, Instrument status is not available. (see Figure 2-9 – Isolation Mode below).


Figure 2-9 – Isolation Mode

- For a stack of BD Phoenix M50 instruments connected to the same AIO PC, each instrument can be in Isolation Mode independent of the other instrument (see Figure 2-10 – Isolation Mode for One Instrument).


Figure 2-10 – Isolation Mode for One Instrument

- The AIO PC handles the transition of each instrument into and out of Isolation Mode independently.
- In Isolation Mode, when you open the instrument door, no station status indicators are lit. Routine workflow is not supported in Isolation Mode.
- The instrument and AIO PC both return to normal operations when communication between the two is reestablished. During the transition, test cycle data collected by the instrument while in Isolation Mode is transferred to the AIO PC and processed. Panel results are evaluated at this time for all panels that are still in the instrument when recovering from Isolation Mode.

Instrument Alert Indicator during Isolation Mode		
LED Appearance	State	Meaning
Blinking	Instrument not connected	The instrument is not connected or it is not communicating with the AIO PC

Isolation Mode Operation

- Routine workflow is not supported in Isolation Mode

Isolation Mode Troubleshooting

Isolation Mode can be caused by the following conditions:

- AIO PC malfunction
- AIO PC power or communication (USB) cable disconnected
- BD Phoenix M50 instrument application on the AIO PC has stopped working

To return to normal operating mode, check USB and power cables and reconnect if needed. Rebooting the AIO PC and/or instrument may also be necessary. If these actions do not correct the problem, contact your local BD representative.

2.8 LIS Operations

2.8.1 General

The LIS communications feature enables the BD Phoenix M50 instrument to exchange information with a compatible Laboratory Information System (LIS). LIS communications can be configured to exchange Order records and Results records at a variety of times.

Order records can be downloaded from the LIS system to the BD Phoenix M50 instrument. These Order records can include the information listed in the second bullet below. If all this information is sent from the LIS system, then the panel is automatically logged in just as if the login were done at the instrument. If the panel sequence number is omitted, the Order record can be associated to a specific panel manually in Panel Login. Order records can be configured to be uploaded to the LIS system also.

2.8.2 Results Upload Records

Results records are uploaded from the BD Phoenix M50 instrument to the LIS system. These records consist of:

- Header record (Delimiter fields, sender name, version number, message date/time)
- Order record (Accession number, Isolate number, Organism, Test ID, Sequence number, Priority, Report type)

Note that if an Organism code is included in the order record, for either a new or existing panel, the ID side of the panel is disabled.

- Comment record (contains Special Messages and/or BDxpert Rules)

2.8.3 Results Record

- Panel sequence number
- Instrument number
- Instrument type
- Instrument location (Station)
- Time to result for identification or MIC produced
- Test start time
- Test end time
- Test status (ongoing, complete, partial complete, complete with needs attention reason, complete with all ignored needs attention reasons, complete QC Pass, complete QC Review, pending, or rapid complete)
- Result type
- Antimicrobial code
- MIC value
- S / I / R / No Interp / Error value
- Resistance marker
- ID or Final ID
- Results status (finalized or unfinalized)

The LIS Communications function is based on the American Society of Testing and Materials (ASTM) LIS Communications Standards (1381 and 1394), and is compatible with a number of popular LIS systems. For specific information on which LIS systems are compatible, contact your local BD representative. The LIS Vendor Interface Specification (available upon request) provides complete details on the BD Phoenix implementation of LIS Communications.

NOTE

If the BD Phoenix M50 instrument is connected to a LIS, it cannot be connected to the BD EpiCenter system. However, if the BD Phoenix M50 instrument is connected to BD EpiCenter, then a LIS connection can be established via BD EpiCenter.

2.8.4 Important Concepts

There are several concepts with which the user should be familiar in order to understand the information presented here, and to ensure that the LIS communications feature is set up properly for your laboratory. These concepts are discussed below.

LIS communication is able to send Results records from the instrument to the LIS (**upload**), order records when panels are placed into the instrument, and queries to (and from) the LIS for Order records. Results upload can be configured to include or exclude Interpretation (SIR) results. The instrument can be set up to upload Results records only when the LIS requests them (**solicited upload**); or at one of the following unsolicited upload times: when panels are finalized; when panels complete testing or when complete panel records change; when ID or AST results are determined; or at a fixed time. QC panels and orphan panels are uploaded only when solicited by the LIS. If the instrument is configured for unsolicited upload, it still responds to requests from the LIS for results (solicitations or queries). If the instrument is configured for unsolicited uploads, the LIS must always be ready to receive data from the instrument.

By downloading barcode (panel) sequence numbers, the instrument can automatically log in panels when records are received. If panel sequence numbers are omitted from the Order record, the panel can still be associated to the Order record in the Panel Login screen.

If the panel sequence numbers are not included in downloaded records, the records that are downloaded are stored (up to 200 Order records) for manual association to panels. This process is described in Section 3.3.2 Inserting Panels in the Instrument.

Organism Configuration and Antimicrobial Configuration enable the user to enter the specific codes required by your LIS system for the organisms and antimicrobials uploaded in Results records. Refer to **Options** in Section 3.10.2 System Configuration Sub-Tab for additional information.

LIS configuration settings are independent of critical panel configuration settings. (e.g., if LIS configuration is set to send results only when the panel is complete, the results are not uploaded if the panel is critical and rapid reporting configuration is set for notification on ID results or partial results.)

2.8.5 Routine System Operation

Operation of the BD Phoenix M50 instrument with an LIS interface differs very little from routine operation of the system. The major difference is the ability to enter panel/accession data via the LIS into the system. With LIS communications, patient information can be logged in at the LIS and transferred to the BD Phoenix M50 instrument. Consult your LIS manufacturer's operation manual for complete instructions on data entry and downloading records.

LIS systems operate either in real-time mode, where the system automatically downloads each Order as it is logged in, or in batch mode where multiple Orders in a group are logged in and downloaded by the user. After patient records are logged in, download them to the BD Phoenix M50 instrument.

Any data sent to the instrument that does not directly correlate to one of the fields defined as the Order record is ignored by the system. Any information sent from the LIS for a Finalized panel is rejected.

After the Order records have been downloaded, and the panels have been attached to those records (if necessary), routine system operation does not differ in any way. For example, the user can continue to perform tasks such as: loading the instrument, printing reports, monitoring the system for complete panels, and performing maintenance; however, it is recommended that the user be especially alert and quickly respond to any system or activity alerts that occur.

3 – User Interface Operations

3.1 Using the Instrument Interface

The touchscreen presents all of the information needed to monitor instrument status, log in panels, set up the instrument, print reports, and perform routine instrument maintenance. For more information on how to use the touchscreen refer to Section 2.3.5 Touchscreen Keyboard. Operations are presented in the form of tabs which, when selected, produce functional, interactive screens, as well as icons that graphically represent the information (e.g., a thermometer indicates the current temperature). The screen's application header presents instrument status information that is updated every few seconds. The middle region of the screen initially presents station status information. Display regions are discussed in greater detail in Section 2 – Controls, Icons, Indicators, and Installation.

The information below details the tab functions (left to right) which appear on the Status screen's application footer.

3.2 Status Tab

The Status tab will show one or two instrument status' that report the following information:

- Removable: number of panels ready for removal
- Ongoing: number of panels in testing progress
- Empty: number of empty panels
- Blocked: available stations that cannot be used due to an error in the system.
- Time: time remaining for test cycle
- Temperature
- Connection status

3.3 Panel Login Tab

Panel Login enables the log in of panel demographic information to be tested. Depending on the type of panel being used, not all fields listed below may appear on the screen. For information on logging in QC panels, refer to Section 4.3.1 Quality Control. More than one of the same panel type may be logged in for an accession number, but the system displays a notification for duplication with an activity error message.

Panel Login Fields:

Accession Number

Type in or scan an accession number, up to 20 alphanumeric characters. If the accession barcode was scanned to access Panel Login, the Accession # is completed automatically. Spaces at the beginning or end of the number are ignored, but spaces within the number are valid. The following characters CANNOT be used in accession numbers: * ? [] ! # |.

Sequence Number

Type in or scan the panel's sequence number. (If the Sequence Number was scanned to access Panel Login, this field is completed automatically.) The sequence number contains digits that identify the panel type. Valid sequence numbers are 12 digits long and must conform to BD panel sequence number specifications. Enter a panel sequence number to save a record. If an existing sequence number is entered, the system displays the Panel Results screen.

Isolate Number

Defaults to isolate number 1. Valid isolate numbers are 1 to 20. If an accession number has been entered, an isolate number must also be entered. If the isolate number is changed from 1 to another value, an accession number must be entered.

Critical Check Box

Defaults to disabled. The system can be configured (Section 3.10.8 Rapid Reporting Sub-Tab) to provide a notification of critical panel results (ID results obtained, partial panel results, or complete panel results) by sounding an audible alarm and/or printing a Lab Report automatically (or neither).

ID Check Box

Defaults to enabled if a combination panel sequence number is scanned. This field appears only when a combination panel is used. This field is automatically checked (enabled) when a QC panel is logged in.

AST Check Box

Defaults to checked if a combination panel sequence number is scanned. This check box appears only when a combination panel is used. This field is automatically enabled when a QC panel is logged in.

Organism ID/Test Strain

Use the drop down list to select the organism ID. Organisms are listed in alphabetical order. Enter the first few characters of the organism name to jump to that portion of the list quickly.

This field appears only for AST-only panels or when using the AST side of a Combination panel.

For QC panels, this field is named Test Strain, and lists only the strains of organisms predefined in the database, sorted by strain number (alphanumerically). A test strain must be entered to save a QC panel.

Tech ID

This field appears for QC panels. Enter the identification for the technologist performing the QC test. Up to 3 alphanumeric characters are accepted. A tech ID must be entered to save a QC panel.

Panel Lot Number

This field appears for QC panels. Type in or scan the panel's lot number. Lot numbers must be 7 digits. A lot number must be entered to save a QC panel. (When the QC Lot Support feature is enabled, this field is completed automatically when the Sequence number barcode is scanned.)

Expiration Date (Panel Lot)

This field appears for QC panels. An expiration date must be entered to save a QC panel. (When the QC Lot Support feature is enabled, this field is completed automatically when the Sequence number barcode is scanned.)

ID Broth Lot Number

This field appears for QC panels. Type in or scan the broth lot number. Lot numbers can be up to 7 characters.

Expiration Date (ID Broth)

This field appears for QC panels. An expiration date must be entered to save a QC panel.

AST Broth Lot Number

This field appears for QC panels. Type in or scan the broth lot number. Lot numbers can be up to 7 characters.

Expiration Date (AST Broth)

This field appears for QC panels. An expiration date must be entered to save a QC panel.

Indicator Lot Number

This field appears for QC panels. Type in or scan the indicator lot number. Lot numbers can be up to 7 characters.

Expiration Date (Indicator Lot)

This field appears for QC panels. An expiration date must be entered to save a QC panel.

3.3.1 Logging in Panels

The inoculum density of the panel is set in Configuration. The density setting cannot be changed during Panel Login. The only way to use a different inoculum density is by blackening well A-17 as described in Section 4.3 Preparing Panels.

In order to insure optimal system performance for Yeast ID panels only, the correct media type must be provided via a drop-down box or by using the default media setting. The media type selection only applies to Yeast ID panels and is not displayed for other panel types. Log the panel into the instrument as follows:

- 1 Select the **Panel Login** tab.
- 2 Select **Clinical**.
- 3 If the receipt of a special notification is desired (audible alarm and/or automatic printing of a Lab Report) when panel results are obtained (ID only, partial, or complete), select **Critical**. More information about critical panels is provided in Section 3.10.8 Rapid Reporting Sub-Tab under **Critical**.
- 4 In the **Accession Number** field, type in or scan an accession number.
- 5 In the **Sequence Number** field, type in or scan the panel's sequence number.
- 6 The **Isolate Number** field defaults to isolate number 1. Type in the isolate number, or touch the +/- to increase/decrease the number. Valid isolate numbers are 1 to 20. Enter an isolate number if an accession number has been entered.
- 7 For Yeast ID panels, a media type must be specified in the **Media** field. If a media type has not been specified, a workflow error is generated when the panel is attempted to be saved. If a

Yeast ID panel is not logged in before placing it in the instrument for testing, the panel aborts after the first reading because no media has been specified. A default media type can be configured (see **OPTIONS** in Section 2 – Controls, Icons, Indicators, and Installation) which appears when a Yeast ID panel sequence number is scanned during login. To select a different media, use the drop down box listing all media types (abbreviations), sorted alphabetically. (When the media type is selected, the full name appears at the top right of the screen.) Select the media type.

- 8 If either the ID or AST portion of a combination panel is only being used, disable the part of the panel you are not using.
- 9 If you disable ID on a combination panel, or if an AST-only or BD Phoenix Emerge panel is not being used, the Organism ID field appears. If the system is not performing the organism identification, an organism ID for SIR interpretation must be provided. (If an AST panel is being tested and an organism ID is not entered, the panel will go to **Needs Attention** when the instrument completes reading. An organism ID will have to be provided in order for the BDxpert system to interpret MIC results.) Highlight the desired organism from the drop-down box. Organisms are listed in alphabetical order. Enter the first few characters of the organism name to move up to that portion of the list quickly. Select the desired organism. The desired organism can also be scanned from the barcode list of organisms found in the Quick Reference Guide.
- 10 Select **Save** to save the information.
- 11 Place the panel in the instrument (refer to Section 3.3.2 Inserting Panels in the Instrument). The following functions can be performed from Panel Login:

Save – saves the information displayed

Repeat Data – enters the last accession number and media types for panels, or media type and lot information for QC panels as follows: if QC Lot Support is disabled (Panel Lot plus Expiration, ID Broth Lot plus Expiration, AST Broth Lot plus Expiration, Indicator Lot plus Expiration), if QC Lot Support is enabled (ID Broth Lot plus Expiration, AST Broth Lot plus Expiration, Indicator Lot plus Expiration)

Cancel – clears the displayed record from the screen

3.3.2 Inserting Panels in the Instrument

Panels are inserted with the reaction (and panel sequence barcode label) sides facing the interior of the instrument.

- 1 Select Panel In (see Figure 3-1 – Inserting Panels)
- 2 When the blue light is blinking, open the instrument door. An audible tone will sound and the unlock icon will be visible.

WARNING

THE INSTRUMENT DOOR IS ELECTROMECHANICALLY LATCHED AND IS CONTROLLED BY THE INSTRUMENT SOFTWARE. NEVER ATTEMPT TO DEFEAT THE DOOR LATCHING MECHANISM, OR TO OPEN THE DOOR WHEN THE UNLOCKED ICON IS NOT DISPLAYED. SERIOUS INJURY CAN BE CAUSED BY THE ROTATING CAROUSEL. IF THE CAROUSEL IS NOT COMPLETELY STOPPED WHEN THE DOOR IS OPENED, IMMEDIATELY CONTACT BD FOR SERVICE. NEVER ATTEMPT TO ROTATE THE CAROUSEL MANUALLY OR SERIOUS INJURY MAY RESULT.

- 3 Select a panel holder where there is no panel in place and no LEDs are illuminated. Place the bottom part of the panel in the panel holder.
- 4 Press downward.
- 5 Pivot the top of the panel back into the panel holder.
- 6 Allow the panel to move upward into place.
- 7 Close the instrument door. If more panels need to be inserted than there are available holders in the current section, wait for a moment for the carousel to rotate to provide additional available holders and repeat **Steps 2 through 7**.
- 8 The system performs an inventory scan to locate any newly inserted panels and reads the bar-codes of these panels.

NOTE

Do not snap the panel back into the holder. This may result in splashing of the inoculum which may cause inaccurate results.


Figure 3-1 – Inserting Panels

3.3.3 Unloading and Discarding Panels

When panel testing is completed, panels should be removed from the instrument and discarded.

WARNING

THE INSTRUMENT DOOR IS ELECTROMAGNETICALLY LATCHED AND IS CONTROLLED BY THE INSTRUMENT SOFTWARE. NEVER ATTEMPT TO DEFEAT THE DOOR LATCHING MECHANISM, OR TO OPEN THE DOOR WHEN THE "UNLOCKED" ICON IS NOT DISPLAYED. SERIOUS INJURY CAN BE CAUSED BY THE ROTATING CAROUSEL. IF THE CAROUSEL IS NOT COMPLETELY STOPPED WHEN THE DOOR IS OPENED, IMMEDIATELY CONTACT BD FOR SERVICE. NEVER ATTEMPT TO ROTATE THE CAROUSEL MANUALLY OR SERIOUS INJURY MAY RESULT.

To remove panels:

- 1 Select **PANELS OUT**.
- 2 When the blue Panel In/Out indicator blinks, open the instrument door.
- 3 All panels that are ready to be removed are indicated by a solid green LED indicator.
- 4 Remove the panels by pushing the panel down, pivoting the top of the panel outward, and pulling the panel out of the panel holder.

-
- 5 If there are completed panels that are not in the accessible stations, close the door and allow the instrument to reposition the carousel to provide access to those panels. Open the door and continue removing completed panels.
 - 6 Discard the panels in a biohazard container.

3.4 Results Tab

Results enables the review and modification of panel test results. Results can be used for the following functions:

- ◆ display a panel whose data is stored in the BD Phoenix database
- ◆ modify the information for a panel in the BD Phoenix database
- ◆ mark a panel as critical
- ◆ print a Lab Report on any panel that can be recalled or displayed on the screen
- ◆ locate a panel resident in the instrument
- ◆ delete panel information from the BD Phoenix database
- ◆ answer or display any triggered BDXpert Rules
- ◆ display any special messages
- ◆ finalize a panel

Results can be accessed in several ways:

- scanning a known panel sequence number or accession number with the external scanner while the Status screen is displayed
- selecting **Results** from the Status screen (see Figure 3-2 – Results Screen)
- scanning or entering a known panel sequence number while the Panel Login screen is displayed
- selecting **Delete** or **AST Results** for a panel in the Needs Attention tab when the Delete option is active
- selecting **AST Results** from the Batch Finalization screen
- selecting **Results** from Inventory (Section 3.7 Inventory Tab) or Panel Lot Definition (Section 3.10.9 Panel Lot Definition Sub-Tab).


Figure 3-2 – Results Screen

Typical Panel Results are shown in the figure below. (Figure 3-3 – Panel Results Screen)


Figure 3-3 – Panel Results Screen

Actions available from all screens are:

- **Locate Panel** – causes the instrument to locate and indicate the current panel and unlock the door
- **Save** – saves any changes
- **Print** – prints a copy of the Lab Report for the current panel
- **Delete** – deletes the panel results
- **Cancel** – returns user to initial Search for Panel Results screen

Panel Results (clinical panels) are retained for 31 days (possibly longer depending on number of QC panels tested). QC panel results are retained for at least six months.

The type of panel is shown in the title area of the screen. The inoculum density used for the identification (if applicable) is shown in the results header area after the first test cycle completes.

BDXpert Rules

When BDXpert Rule is selected, the following are displayed (note that all fields may not be displayed depending on the type of panel and whether or not the instrument is connected to a BD EpiCenter system):

- Rule Number
- Status
- Resistance Marker
- Name
- Rule Description

ID and Biochemicals

When ID and Biochemicals is selected, the following appears:

- Instrument Organism ID which may consist of Instrument Organism ID, Confidence Value, and Up to 5 Supplemental Tests if the panel was set up using low inoculum
- Biochemical Results which consist of: Biochemical Abbreviation, Actual Result (+, -, ?, or X) or Expected Result (+, -, V, or blank) shown in Figure 3-4 – ID and Biochemical Results
- Special Messages (icon represented in Figure 3-4 – ID and Biochemical Results)
- Needs Attention (icon represented in Figure 3-4 – ID and Biochemical Results)


Figure 3-4 – ID and Biochemical Results

If a Needs Attention (see Section 3.6 Needs Attention Tab) exists for the panel, it appears in the results header area. If a special message (see Section 3.4.8 Special Messages Sub-Tab) exists for the panel, an icon is shown in the results header area. If a Resistance Marker has triggered for the panel, the Resistance Marker icon (see Section 2.5 Software Icons) is shown in the results header area.

If a panel is recalled by accession number only, and there is more than one panel attached to the accession, the Results List screen appears. From this screen, select the specific panel to review/modify.

Note that depending on the type of panel being used, as well as other circumstances, not all fields listed below may appear on the screen.

AST or Combination panel types (with at least the AST side enabled) each contain a set of antimicrobials. The instrument reports a result for each antimicrobial on the panel. Each antimicrobial reports an individual MIC value. Once an antimicrobial has a MIC value, the instrument calculates an Instrument susceptibility (SIR) value for each MIC value that has been determined. (The instrument requires the panel to have an Organism ID defined to interpret MIC values into Instrument SIR values. Also note that SIR values are not calculated by the instrument if it is attached to an EpiCenter Data Management System.)

If Rapid Completion is enabled, the instrument is able to provide BDXpert AST results (SIR) before determining actual MIC values. The instrument MIC values are provided as soon as they can be accurately determined. Within a test panel, some MIC values may be available earlier than others. The rapid completion feature can be used to predict resistance for uncompleted antibiotics using the ID alone (intrinsic resistance), or ID with completed MICs of related antibiotics, or resistance marker tests (BL, ESBL). The BDXpert system is used to make these predictions. This can be useful in situations where, for example, the results for drugs that have not yet received MICs would be of no clinical value based on the other results that are already available. Antimicrobials with Rapid Complete BDXpert interpretations are indicated by a C in the MIC column on Results screens and Lab reports.

When both a MIC value and Instrument SIR value have been determined for an antimicrobial, the instrument executes the BDXpert Rules (providing the BDXpert System is enabled). The instrument

reports a value in the BDxpert SIR field if an enabled BDxpert Rule triggered and the reported BDxpert SIR value is different from the value in the Instrument SIR field.

Different results appear depending on which tab is selected. The header information remains the same no matter which tab is selected:

Panel Results Fields:

NOTE

Modifying the Accession Number or Isolate Number may invoke Auto Association which can change results.

Accession Number

Type in or scan the accession number to recall. If only an accession number is entered, and there is more than one panel attached to the accession, the Results List screen appears.

This field can be modified for unfinalized panels. Enter up to 20 alphanumeric characters for the accession, excluding * ? [] ! # |. Modifying an accession number does not affect the accession number of any related panels.

Sequence Number

This field cannot be edited.

Isolate Number

Valid isolate numbers are 1 to 20. This field can be modified for unfinalized panels, however an existing isolate number cannot be changed to a blank number.

Media

This field is editable only before the first test cycle completes.

Shows the media type selected during Panel Login for Yeast ID panels only. The following values can appear in this field: blank for Unspecified Media Type; INVLD (Invalid Media Type); SABDX (Sabouraud Dextrose Agar); TSASB (BD Trypticase™ Soy Agar w/5% SB); COLSB (Columbia Agar w/5% SB); CHOC (Chocolate II Agar); SABEM (Sabouraud Dextrose Emmons); SABHI (Sab Brain Heart Inf Ag Deep)

Status

Read-only field shows the panel's testing status: Pending, Ongoing, Complete, or Rapid Complete (if enabled). Rapid Complete panels have not finished testing and show BDxpert results. If a Rapid Complete panel is removed from the instrument, it then becomes Complete. If it is left in the instrument, Rapid Complete MICs are replaced by actual MIC results as they are determined.

Final ID

Final ID can be completed automatically by the system from a single instrument-based ID or by the user selecting among the instrument-based IDs. For ID panels (or ID-only portions of Combination panels) organism ID cannot be manually selected if there is no instrument-based ID listed.

Note that if the **Modify Related Panels** check box is checked, the ID in the database of all related

panels (panels with the same accession and isolate number) is set with the same organism ID through the Auto Association function (see Section 4.4 Automatic Association of Panels).

Organisms are listed in alphabetical order. Enter the first few characters of the organism name to find that portion of the list quickly.

When the field is modified, the system re-evaluates each antimicrobial's instrument SIR value, as well as re-evaluating BDxpert information.

The field cannot be modified if the panel is finalized.

For QC panels, this field is named **Test Strain**, and lists only the ATCC strains of organisms predefined in the database, sorted by strain number.

Panel Lot Number

Read-only field shows the panel lot number. This field shows only when the QC Lot Support feature is enabled.

Expiration Date

Read-only field shows the date and time that the panel lot expires. This field shows only for QC panels and only when the QC Lot Support feature is enabled.

Modify Final ID of Related Panels Check Box

This field appears when the Final ID field is modified. Defaults to checked when the Final ID is modified. If the ID for unfinalized related panels is NOT to be modified when the currently displayed panel is modified, disable this field as it has no affect on related QC panels. This field does not appear if the EpiCenter system is attached and communicating.

Test Start

Read-only field shows the date and time that panel testing was started.

Test End

Read-only field shows the date and time that panel testing was completed.

Location

Read-only field shows the location of the panel, in the form Tst, where T is the Tier, and st is the station number.

QC Status

This field shows the status of a QC panel. The field is blank until the panel status becomes Complete. Statuses are initially PASS or REVIEW. Review indicates that the panel has not passed. Check any panels with Review status and determine why the panel did not pass.

The status is REVIEW if any of the following conditions occur:

- ◆ QC strain was identified incorrectly
- ◆ The test on a QC panel is aborted
- ◆ At least one of the antimicrobial results fail

From a status of REVIEW, the final status may be set to REPEAT if it is determined that the panel failed due to error in preparation or handling of the panel. If it is not determined that a panel preparation/handling error was made, the final status should be set to FAIL. Selecting REPEAT or FAIL clears the Review QC Results Needs Attention condition.

Finalized

Check the box to finalize the panel. The **Finalized** field is displayed as a read-only field when the BD Phoenix system is connected and communicating with the BD EpiCenter Data Management Center. In this case, all Finalization is done at the BD EpiCenter system. **Finalized** is not displayed for QC panels.

If a panel with Rapid Complete status is finalized, MIC results processing stops, and the panel status becomes Complete. Any drugs that did not complete testing maintain their current MIC value (e.g., C in MIC column remains C, X remains X, etc.).

Critical

Enable this field to mark the panel as critical. An audible alarm and/or an automatically printed Lab Report (or neither) notification can be configured (Rapid Reporting configuration (Section 3.10.8 Rapid Reporting Sub-Tab) for critical panel results (ID results have been obtained, partial panel results, or complete panel results).

Previously marked critical panels can be disabled.

This field does not appear for QC panels.

3.4.1 Adding/Modifying ID Results

The system allows an organism ID to be entered manually, or the organism identified by the instrument can be overridden. In addition, in some cases, the system will not be able to make a single identification determination based on panel results. In these cases, two or three organisms may appear in the Instrument ID field. When more than one organism appears as the Instrument ID, the system does NOT automatically enter an identification in the Final ID field. The user must select the desired organism. The actual organism may be determined either through supplemental tests, which are recommended in the Instrument ID window and/or the Special Messages screen (accessible via **Special Messages**), or through other tests performed. To add/modify the Final ID:

- 1 Select Final ID field in the Results header.
- 2 Highlight the desired organism.
- 3 Press the Save button to save.

The following icon appears:

3.4.2 Modifying AST Results

There may be times when the Final SIR results for a panel need modification. Note that Final SIR results cannot be modified if there are manual BDXpert rules pending. The pending rules must first be accepted or rejected, which allows the system to perform its final results processing. After the final processing is complete, the SIR results can be modified manually if desired.

To modify the Final SIR results:

- 1 Select **AST Results** from Results
- 2 Select the **FINAL SIR** field for the desired antimicrobial. The following selections are available:
 - S(usceptible)
 - I(ntermediate)
 - R(esistant)
 - X = Invalid, cannot interpret
(see Table 3-1 – MIC/SIR Values and Causes in Section 2 – Controls, Icons, Indicators, and Installation)
 - N(not Susceptible)*
 - Blank (indicates ID is required)

* N indicates that the antimicrobial/organism does not have an upper breakpoint. When there is no upper breakpoint there are no criteria for calling an organism intermediate or resistant. This often occurs when there are no known resistant strains of an organism. In this case, if the MIC is below the lower breakpoint the SIR results can be reported as susceptible but if the MIC is above the lower breakpoint the only result that can be reported is N or not susceptible.

3.4.3 Modify Panel Usage Sub-Tab

Results will display Modify Panel Usage when the currently recalled panel:

- is an ID/AST Combination panel type and both sides of the panel are enabled
- status is not complete
- has no instrument organisms determined for the panel
- does not have any AST Complete set on its AST side

The user can disable the ID or the AST side of the panel when no other changes have been made to the panel on the Results screen. When the user disables the ID or AST side of the panel, no other tab will be accessible. If the user attempts to leave the tab, a message will be displayed instructing the user to save or cancel the panel usage change.

To modify panel usage:

- 1 Select **Results**.
- 2 Enter or scan the panel barcode sequence number of the panel whose usage you wish to modify. The system automatically completes the Accession # and Isolate # fields, which are read-only.
- 3 Select **Modify Panel Usage**.
- 4 The following screen appears:

Figure 3-5 – Modify Panel Usage Screen

- 5 Select **ID** or **AST** for the side of the panel to be disabled. If both are selected then workflow alert code W305 will be displayed. The screen data will be maintained and the save attempt will be stopped.
- 6 Select **Save** to save the panel modification.

Note the following conditions for panel usage modification:

- No information for related panels is modified when panel usage is changed.
- If there is no user-entered Final ID and the ID side is disabled, the Instrument ID, Biochemical Results, Confidence Values, SIR Values, and ID Special Messages are removed from the record.
- If there IS a user-entered Final ID and the ID side is disabled, this ID is retained, as are SIR Values. However, Instrument ID, Biochemical Results, Confidence Values, and ID Special Messages are removed from the record.

3.4.4 AST Results Sub-Tab

Antimicrobial

Read-only field showing the abbreviation and name for the antimicrobial.

MIC

The Minimum Inhibitory Concentration value determined by the instrument. The following values may also appear in this field:

- > growth occurred for all concentrations of the antimicrobial
- ≤ no growth occurred in all of the concentrations of the antimicrobial
- ? MIC determination is pending (SIR values remain blank)

- X MIC value cannot be produced; or Final ID is not claimed in the Taxa listing (see Section 7.6 Taxa for ID/AST Determination) for AST testing; or the panel's drug dilution series does not cover the BD Phoenix reportable MIC range, or for additional causes, see Table 3-1 – MIC/SIR Values and Causes). Check the Special Messages displayed on the screen (Section 3.4.8 Special Messages Sub-Tab) for an explanation.
- C Rapid Complete (MIC is pending, BDxpert SIR is based on ID and completed drug and/or resistance marker results). These values are replaced by actual MIC values as they are determined.

For QC Panels, the following values may appear in this field:

Actual MIC value (number, ?, or X for error)

Expected MIC range (if defined)

Status

- P the actual MIC value is within the expected MIC range
- F the actual MIC value is not within the expected MIC range
- R repeat: the actual MIC value is X
- blank no range is defined for the antimicrobial/organism combination

I(nstrument SIR)

This field is the instrument interpretation for the MIC based on the breakpoints currently running in the system (selected in Section 3.10.10 Instrument Sub-Tab).

The SIR value is blank for antimicrobials that require an ID to perform the SIR interpretation. In this case, the **BDxpert SIR** field (if present) and the **Final SIR** field are also blank.

The following represent the interpretation values:

- S Susceptible
- I Intermediate
- R Resistant
- N Not susceptible*
- X Cannot produce interpretation

Blank No SIR (missing or invalid Final ID; MIC = ?, C, or X; unclaimed organism for Final ID)

* N indicates that the antimicrobial/organism does not have an upper breakpoint. When there is no upper breakpoint there are no criteria for calling an organism intermediate or resistant. This often occurs when there is an absence or rare occurrence of resistant strains of an organism. In this case, if the MIC is below the lower breakpoint the SIR results can be reported as susceptible but if the MIC is above the lower breakpoint the only result that can be reported is N or not susceptible.

I(nstrument SIR) does not appear if the BD EpiCenter system is attached and communicating. If communications with the BD EpiCenter system is lost, the field automatically appears.

(BD)E(Xpert SIR)

Possible field values are: Blank, S, I, R, N, and X. These values have the same meanings as I(nstrument SIR) above.

The BD(Xpert SIR) field is not shown when the BD(Xpert System) is disabled. It is also not shown for QC, ID, or ID/AST panels with AST disabled. This field contains the results calculated by the BD(Xpert System) based on the execution of all enabled BD(Xpert Rules). Values appear in this field only if a BD(Xpert rule) triggered and caused the BD(Xpert SIR) value to differ from the instrument interpretation.

Rule

Read-only field that shows the numeric designation of the highlighted rule.

F(inal SIR)

Field values are the same as I(nstrument SIR) above. The Final SIR can be one of the following (from highest to lowest priority): a user-entered SIR value; the BD(Xpert SIR) value; the instrument SIR value; or blank.

Final SIR values are produced when a MIC value or error is determined, AND there are no pending manual BD(Xpert rules (if the BD(Xpert System) is enabled).

When the BD EpiCenter System is not connected and the BD(Xpert System) is disabled, the Final SIR field is not displayed for an antimicrobial until the MIC value is something other than ?.

The following table provides explanations of different combinations of blank and X results. Detailed explanations for actual results are provided as Special Messages (Section 3.4.8 Special Messages Sub-Tab).

Antimicrobial	MIC	Instrument SIR	BD(Xpert SIR)	Final SIR	Possible Cause Examples
Any Drug	?	[Blank]	[Blank]	[Blank]	Results pending.
Any Drug	BD Phoenix MIC Result	S, I, R	[Blank]	S, I, R	Instrument SIR = Final SIR (No BD(Xpert rule SIR).
Any Drug	BD Phoenix MIC Result	S or I	R	R	Final SIR = BD(Xpert SIR) if rule accepted.
Any Drug	BD Phoenix MIC Result	[Blank]	[Blank]	[Blank]	No breakpoints for this drug/organism combination within the chosen standard (CLSI, SFM, EUCAST, Custom).
Any Drug	BD Phoenix MIC Result	X	[Blank]	X	MIC value is outside the breakpoints for the selected standard. Example: Panel drug range = 1–16 µg/mL, susceptible breakpoint = 0.5 µg/mL
Any Drug	BD Phoenix MIC Result	X	[Blank]	X	SIR is suppressed by a BD(Xpert rule). User must provide Final SIR based on manual interpretation or additional testing.

Antimicrobial	MIC	Instrument SIR	BDXpert SIR	Final SIR	Possible Cause Examples
Any Drug	X	[Blank]	[Blank]	[Blank]	The MIC for this drug and organism combination is not reported by the BD Phoenix system. An alternative method should be used.
All Drugs	X	[Blank]	[Blank]	[Blank]	This species is not included in the BD Phoenix AST taxonomy; perform an alternative method. An excessive amount of indicator was detected in the panel. The AST portion of the panel has been terminated and the isolate should be retested.
Any Drug	X	[Blank]	[Blank] or R	[Blank] or R	The MIC for this antibiotic is not reported (see Special Message). Interpretation based on BDXpert Rule.
Any Drug	C	[Blank]	R	R	Rapid Completion SIR (BDXpert SIR based on ID and/or another completed drug and/or resistance marker result).

Table 3-1 – MIC/SIR Values and Causes

3.4.5 BDXpert Rules Sub-Tab

The BDXpert Triggered Rules screen provides a view of the BDXpert system rules that have been triggered for a panel (the panel currently selected in the Results tab). Other views available are: a listing of rules that have been triggered, the text of those rules, the effect that the rules have on Final SIR values, and the ability to accept or reject pending (manual) rules. Additionally, all rules can be re-run.

After all rules have been reviewed and Accepted/Rejected, any changes must be saved.

The BDXpert Rules tab does not appear, and the screen is not available if the BD EpiCenter system is attached and communicating with the BD Phoenix M50 instrument. However, if communications with the BD EpiCenter system is lost, and the BDXpert system is reactivated in Instrument Configuration (Section 3.10.10 Instrument Sub-Tab), the tab reappears and the screen can be accessed. When communications with the BD EpiCenter system is restored, BDXpert rules interpretations are once again performed at the BD EpiCenter system.

Rules Field

Rule #

Read-only field that shows the numeric designation of the highlighted rule.

Status

Shows the status of the rule. The initial status of Automatic (a rule that executes automatically without user intervention) or Manual (a rule that must be manually accepted or rejected) is set in the BDXpert Rule Configuration screen (Section 3.10.7 BDXpert Rules Sub-Tab). Statuses are:

Automatic – rule is enabled and set to Automatic

Pending – rule is enabled and set to Manual; Manual rules must be Accepted or Rejected; only the first Manual rule shows as Pending

Accepted – rule is enabled and set to Manual and has been Accepted by the user

Rejected – rule is enabled and set to Manual and has been Rejected by the user

Pending rules can be accepted or rejected via **Accept** or **Cancel**. Once a rule is accepted or rejected, the status can only be changed by re-running the rules.

3.4.6 Lot Information Sub-Tab

This sub-tab provides a listing of disposables used in the set up of panels. To access this sub-tab, follow the steps below:

- 1 Enable QC Lot Support (see Section 3.10.2 System Configuration Sub-Tab) in System Configuration.
- 2 Go to Panel Lot Definition (see Section 3.10.9 Panel Lot Definition Sub-Tab) and scan in the barcodes.
- 3 Go to Results and enter an Accession Number.
- 4 Select Search.
- 5 Lot Information Sub-tab appears (see Figure 3-6 – Lot Information Sub-Tab Screen).

The screenshot shows a software interface for lot information. At the top, there's a header bar with icons and a title bar. Below the header, the main area is divided into several sections. The top section contains fields for 'Accession Number', 'Sequence Number', 'Isolate Number', 'QC Status', 'Test Entry', 'Status', 'Location', 'Inoculum Density', 'Test Start', 'Test End', and 'Tech ID'. The 'Status' field is set to 'PENDING'. Below this, there's a section for 'Lot Information' with fields for 'Panel Lot Number', 'Expiration Date', 'AST Broth Lot Number', 'ID Broth Lot Number', 'Indicator Lot Number', and 'Expiration Date'. At the bottom, there's a toolbar with buttons for 'Search', 'Print', 'Export', and a bottom navigation bar with icons for 'Home', 'Panel Setup', 'Results', 'QC', 'Report', 'Inventory', 'Reports', 'Maintenance', and 'Configuration'.

Figure 3-6 – Lot Information Sub-Tab Screen

3.4.7 Instrument Organism ID Sub-Tab

Read-only field showing the organism identified by the instrument based on biochemical results in an ID or ID/AST panel. Organism IDs listed for a Tie (two possible organisms meet the results) are listed in the Instrument Organism ID field as **organism 1/organism 2**. In some cases the system will not be able to make a single identification determination based on panel results. In these cases, two or three organisms may appear in this field. The user must select the desired organism to enter into the Final ID field.

Confidence

The Confidence value computed by the system is based on the actual biochemical results versus the expected results. The Confidence value is a percentage from 0 to 99.

Supplemental Tests

Supplemental Tests are displayed in the Instrument Organism ID field and/or on the Special Messages screen if there is more than one Organism listed there. Once these tests have been performed, the results of the tests will help to distinguish which organism ID to associate to the AST results. At this point, a single organism ID can be selected from the Final ID field.

Biochemical

Read-only field showing the abbreviation for the biochemical used to determine the ID.

Actual

Read-only field showing the observed biochemical result at the time when the organism ID was determined: + for positive or – for negative. A ? indicates that the biochemical result is pending or that the test was aborted prior to ID results being determined. The letter X indicates an error.

Expected

Read-only field showing the expected biochemical result according to the Instrument Organism ID: + or – for the organism. The letter V indicates that the result can be variable. This field is blank when the actual results are ? (for non-QC panels), or until an Organism ID has been determined, or when more than one organism is listed in the Instrument Organism ID field.

Resistance Marker Field

Resistance markers are shown in the order in which they trigger. The following information is shown when a Resistance Marker is detected and is displayed on the BDxpert Rules tab:

- The BDxpert rule number that triggered the Resistance Marker
- The Resistance Marker code (abbreviation)
- The Resistance Marker Name
- The BDxpert rule Description

3.4.8 Special Messages Sub-Tab

Special Messages allows information to be accessed about certain panel ID and AST results, as well as some panel readings. These special messages are triggered and are available for viewing regardless of whether the BDxpert System is enabled or disabled. If a recalled panel has an associated **Special Message**, an icon is displayed in the header area of the results screen.

Special Messages on a recalled panel are shown according to hierarchy on the Special Messages window.

The Special Messages screen reflects messages that exist at the time the screen is accessed: it is not updated dynamically with messages that are triggered after the screen has been accessed. To view newly triggered messages, recall the panel again, and re-enter **Special Messages** again.

3.4.9 Needs Attention Sub-Tab

Needs Attention provides a list of all needs attention reasons, listed in priority order as shown in Table 3-2 – Needs Attention Resolutions. When any of these conditions occur, the panel will have a Needs Attention reason code set and the Panel with a Needs Attention reason will be displayed on the Needs Attention screen. To ignore a Needs Attention reason, select the corresponding field next to the condition. For additional information on Needs Attention see Section 3.6 Needs Attention Tab.

3.4.10 Lab Report

The Lab Report can be printed from the Results or Finalization screens and is not available from the Reports tab. It contains all information for a panel sequence number that exists in the BD Phoenix database, including all information in the Results screen, any special messages, BDXpert Rules that triggered, or Needs Attention reasons if they exist.

The QC Lab Report is accessible from the Reports menu and Results screen. It provides similar information to the Lab Report, but prints when a QC panel is being displayed and a report is requested.

The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), and Date and Time Printed, Software version/PUD version.

Body of Report: Top Region: Accession Number, Isolate Number, Sequence Number, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Lab Report is sorted by Accession # then by Isolate # within accessions.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDXpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDXpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report.

ID and AST sections of combination panels print on separate pages of the report.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

The system can be configured (Section 3.10.10 Instrument Sub-Tab) to print an abbreviated lab report. The abbreviated report does not contain the individual biochemical results for ID or ID/AST panels.

3.5 Finalization Tab

Finalization enables the review and finalization of panel test results in a batch mode or individually. Panel results can be finalized for one or all panels meeting the specified conditions, skip to a new panel, or modify and finalize panel results for a single panel. All fields on the Finalization screen are read-only. To modify information for a panel, select **Results**. When the instrument is connected to the BD EpiCenter Data Management System, the Finalization screen is not available.

When **Finalization** is selected, the instrument finds the records that are eligible for finalization. To be eligible, a panel must have a status of Removable, not be a QC panel, and have no unignored Needs Attention reasons. Eligible panels are sorted first by accession number and then by isolate number.

When **Batch Finalize** is selected, the hourglass icon indicates that the finalization operation is in progress. A maximum of 200 panels may be finalized in a given finalization session.

To finalize/batch finalize panels:

- 1 If there are no panels to be finalized, no data appears on the Finalization screen. If there ARE panels to finalize, the following options appear:
- 2 To finalize ALL eligible panels, select **Finalize**.
- 3 To finalize panels one at a time, check the box next to the panel, then select **Finalize**.
- 4 To print a summary report of all panels eligible for finalization, select **Finalized**. This report shows the Accession Number (primary sort), Isolate Number (secondary sort), Test End date and time, Sequence Number, Instrument Number/Station, and Finalized status (* if finalized, blank if not). Select the panel row to go to the **Results** tab.
- 5 If **Results** is accessed to add or modify information before finalizing a panel, be sure to save the modifications, and return to the **Finalization** screen.
- 6 Once an individual (or batch of) panel(s) has been finalized, the **Finalized** button appears on the Finalization screen. This action enables the printing of Lab Reports for all panels that have been finalized during this session (up to 200 maximum).
- 7 Continue to review panel records and finalize until no additional panel records are displayed.

Two types of reports, standard Lab Report(s) and the Finalization Summary Report, may be printed. The Lab Report can only be printed after one or more panels have been finalized. The Summary report may be printed any time.

Finalization fields:

Sequence Number

Read-only field showing the panel's sequence number.

Accession Number

Read-only field showing the panel's accession number.

Isolate Number

Read-only field showing the panel's isolate number.

Media

Read-only field showing the media type selected during Panel Login for Yeast ID panels.

Status

Read-only field showing the panel's status.

Final ID

Read-only field showing the panel's final organism ID.

Test Start

Read-only field showing the panel's start of testing date and time.

Test End

Read-only field showing the panel's end of testing date and time.

Location

Read-only field showing the panel's location.

Finalized

Read-only field showing the panel's finalized status.

Instrument Organism ID

Read-only field showing the instrument identification of the organism.

Confidence

Read-only field showing the confidence value for the Instrument Organism ID.

Biochemical

Read-only field showing the biochemical abbreviation.

Actual

Read-only field showing the actual biochemical result.

Expected

Read-only field showing the expected biochemical result.

Antimicrob(ial)

Read-only field showing the drug name.

MIC

Read-only field showing the Minimum Inhibitory Concentration value.

Interpretation

Read-only field showing the Susceptibility interpretation.

BDXpert SIR

Read-only field showing the final results calculated by the BDXpert System based on the execution of all enabled BDXpert Rules as if they had been executed automatically.

Final SIR

Read-only field showing the Final Susceptibility interpretation. To modify a Final SIR value, select **Results**.

3.6 Needs Attention Tab

Needs Attention displays a list of panels in the instrument's database that have encountered a condition that requires operator attention. These conditions generally represent problems with the panels themselves, or with the information related to the panels. In many cases, the conditions can be corrected by operator action.

When viewing the Needs Attention panel list, the system provides the opportunity to resolve or ignore the condition that caused the panel to be placed in the list. If the panel was placed in the list due to missing or unresolved information (e.g., a tie), the instrument provides the ability to add or edit the information to resolve the condition. If the panel was placed in the list due to a software, panel, or hardware error, the instrument provides the ability to delete the panel to resolve the error condition. Deleting a panel that is still testing causes that panel's protocol to be aborted. Only panels whose Needs Attention reason has not been ignored are shown in the screen.

The Needs Attention screen lists the highest priority reason (see Table 3-2 – Needs Attention Resolutions below). Selecting a panel opens the results screen in the Needs Attention tab, listing additional reasons if they exist. Two reports are available for printing: all Needs Attention reasons, and highest priority Needs Attention reasons.

Condition	Active Operations
Test Aborted	• ignore • locate panel • delete panel
Can Not Identify Barcode	• ignore • locate panel
Can Not Read Panel Wells	• locate panel • delete panel
Panel Lot Expired	• ignore • clear NA
Invalid AST Results	• ignore
Panel Missing	• ignore • delete panel
No Growth on Panel	• ignore • locate panel • delete panel
Panel Lot Undefined	none to address condition
Review QC results	• ignore • locate panel • delete panel • panel results
Missing Accession Number	• locate panel • panel results
Missing Organism ID	• panel results

Condition	Active Operations
Can Not Determine Organism ID	• ignore • panel results
Invalid Organism ID	• ignore • panel results
Organism ID Conflict	• panel results
BDXpert Rule Flagged (if enabled and manual)	• panel results
Pending Too Long	• ignore

Table 3-2 – Needs Attention Resolutions

The Needs Attention screen lists the first 100 panels with a Needs Attention, sorted by reason code (in the same order as the list above) and then by accession number within each reason code. Panels without an accession number are listed first within each Needs Attention reason code. Note that even after panels have been removed from the Needs Attention list, the highest priority reason code in the Panel Results screen can still be viewed.

Needs Attention resolution options

Ignore Check Box – check to Ignore the Needs Attention condition. Ignore does not correct the Needs Attention reason code, but it informs the system that the problem has been acknowledged by the user.

Reason – lists the reason for the Needs Attention.

Locate panel – causes the carousel to rotate to the panel's location and lights the station where the panel resides.

Save – saves any selections that are made.

Print – enables the printing of a Needs Attention List report.

Delete – delete all panel results.

Cancel – removes any selections that have been made.

These options are designed to enable the correction of (whenever possible) the condition causing the panel to need attention. Where it is not possible, other resolutions are provided. **AST Results** displays the Panel Results screen, where the Needs Attention condition may be able to be corrected by adding information to the panel record.

To resolve panels that need attention:

- 1 Access the Needs Attention screen.
- 2 Refer to the chart below for detailed information on the particular Needs Attention reason.
- 3 Press the AST Results tab to access Results (to add or modify information), the Delete button to delete the panel, check Ignore or Locate Panel to find the panel in the instrument.

Condition	Meaning	Possible Cause(s)	Resolution(s)
Test aborted	A condition occurred which caused the panel to be invalid	• Ongoing panel not tested for more than 1 hour • Instrument turned off for more than 1 hour • Instrument door open more than 1 hour • Panel moved to a different tier/instrument • Incubator temperature too high or too low • System software did not execute testing algorithms for more than 1 hour • Media type not specified for Yeast ID panel	• Delete the panel • Repeat testing
Cannot identify barcode	Internal barcode scanner could not read a panel barcode in a station where the instrument could determine that a panel was present	• Barcode label obscured or missing • Unknown panel type was placed into the instrument	Locate the panel in the instrument and examine the barcode • If the barcode is obscured, the panel must be discarded and another inoculated • If the barcode appears to be intact, replace the panel and close the door. After the next inventory, check the Panel Needs Attention screen. If the panel does not appear, the internal scanner can now read the panel
Cannot read panel wells	Internal barcode scanner has read a sequence number in a station but the instrument does not detect that a panel is present in that station	Panel not seated properly	Locate panel, remove it and replace it, as above
Panel lot expired	A panel was logged in (or has a test start date) with a panel lot number that has already expired	Panel being logged in or placed in the instrument is from an expired lot	Discard panel, reinoculate isolate using unexpired panel lot
Invalid AST results	At least one MIC cannot be interpreted. (Excludes QC panels.)	Refer to Table 3-1 – MIC/SIR Values and Causes	Repeat testing of the antibiotic that cannot be interpreted
Panel missing	Internal barcode scanner read a sequence number on an Ongoing panel, but the panel is missing	• Panel removed before the test was completed • Internal scanner failure and the sequence number can no longer be read	If the panel is replaced on the same tier within 1 hour after removal, testing will resume. If the panel is not replaced, testing will abort
No growth on panel	No growth in growth control well. (Excludes ID only and QC panels.)	Instrument did not detect growth in the growth control well of the panel	• Subculture the organism (to insure that it is viable) and inoculate a new panel. Panel has been aborted • From the Needs Attention screen delete or ignore the panel
Panel lot undefined	A panel has been entered whose lot number is undefined (non QC panel)	Panel is from an undefined lot	Define panel lot or ignore
Review QC results	Status of a completed QC panel is "Review"	QC panel which yielded an incorrect ID or incorrect AST result for at least one antibiotic, or has no growth in growth control well	Repeat QC organisms. Check: culture purity, inoculum density
Missing accession number	Panel is missing accession or isolate information. (Orphan panel.)	Failure to enter accession or isolate information	Press the "panel results" key. Add the accession number using the barcode reader or by typing it in
Missing organism ID	Panel has no organism ID. (ID required to determine SIR results. Excludes QC panels.)	• For AST only panel, no ID has been entered • The panel has an unresolved tie or triplet instrument ID and has no related panel with an ID	Select the organism ID. Any BDxpert rules triggered by the given ID will automatically be presented at this point. Select or ignore the rules and save. When completed, exit to see complete test results.

Condition	Meaning	Possible Cause(s)	Resolution(s)
Cannot determine organism ID	Panel has an Instrument ID of "No Identification" or has a related panel with "No Identification" as a final ID	Panel has been in test for 12 hours and the instrument cannot determine the identification	Repeat testing. Check the following: • Culture purity • Inoculum density • Correct panel used? • Organism may not be in the BD Phoenix database
Invalid organism ID	Organism ID is not in BD Phoenix database	Panel received download information of an organism ID that is not in the BD Phoenix database	Use alternate method
Organism ID conflict	Panel completes testing and has at least one related unfinalized panel that contains a different Final ID. Excludes QC panels.	Completed panel has at least one related panel that contains a different ID	Select Results and choose an organism. Selecting the organism may trigger BDExpert rules. If the rules are configured as manual, BDExpert appears
BDExpert Rule flagged	Panel triggered at least one BDExpert rule and the rule is manually enabled in Configuration.	A BDExpert rule needs to be invoked in order to determine AST results. A panel is flagged if BDExpert rules are configured as enabled/manual. (Rules that are configured as enabled/automatic will automatically "trigger" and the panel will not be displayed in Needs Attention.	Each Expert rule is displayed individually in sequence. To accept the rule, select Accept. To reject the rule, select Reject. Use ReRun to delete BDExpert system decisions and start over. When all rules have been displayed and dispositioned, complete AST results show. After all rules have been invoked, select Special Messages (if present) to view Special messages about characteristics of the organism.
Pending too long	Panel has not been scanned (during an inventory count) within 30 minutes of logging in Panel Login.	Panel was logged into the instrument but was not placed in the instrument within 2 reading cycles (approximately 30 minutes).	Repeat testing. Delete the panel

3.7 Inventory Tab

Inventory provides a listing of all panels in the instrument (except temperature reference panels). This listing can be sorted in ascending or descending order by the following fields: sequence number, accession number, results, or needs attention reason. The initial default sort order is by accession/isolate number in ascending order. Subsequently, the list defaults to the last sort criteria and screen configuration (primary/secondary) used. If there are no panels in the instrument **No Data Available** appears on the screen.

The top of the Panel Inventory screen shows the number of panels in the list (if there is more than one panel), the sort field, and whether the sort field is ascending or descending.

The following fields are shown on the Panel Inventory screen (if the information is known):

- A ! (exclamation mark) at the beginning of the row indicates that a critical panel or a panel with a resistance marker has not been acknowledged.
- Critical - checked if panel has been marked as Critical at login.
- Sequence Number (of the panel).
- Accession Number (panels without accession numbers are listed first.)
- Isolate Number.
- Inoculum Density (blank for AST panels; ? for ID panels until first test completes).
- Status (ongoing; complete; rapid if Rapid Completion is enabled and has been triggered for a panel).
- Final (Organism) ID.
- Results (**final** if panel is complete and there are no active Needs Attention conditions) [all MIC results are determined for an AST panel, or all MIC values *and* the organism ID are determined for a Combination panel, or the organism ID is determined for an ID panel]; **partial** if a panel is ongoing or complete but has an active (unignored) Needs Attention condition [at least one MIC value is determined for an AST or Combination panel, or the organism ID is determined for a Combination panel]; **none**, if no MIC values or organism ID is determined for any type of panel.

- Needs Attention - the active Needs Attention icon appears if an unignored Needs Attention reason exists; **blank** when no Needs Attention reasons exist or they have all been ignored).
- Instrument - in which instrument the panel is located.

After the list appears, highlight a panel and access the Results screen to view or edit panel information, perform an instrument locate panel operation, or to print a Lab Report for panels with final or partial results.

3.8 Reports Tab

Numerous instrument reports are available for printing:

- Accession Lab Report – provides information for a specified accession number.
- Antimicrobial Code Report – displays all antimicrobials, BD and LIS codes for all antimicrobials, that exist in the BD Phoenix database from all panel configurations.
- BDxpert Rule Set Database Report – lists each BDxpert rule number and the text describing the rule, whether each rule is enabled/disabled and whether each rule shall trigger automatically/manually in the system.
- Completed Lab Report – contains information for all panels whose status became complete during the selected time period (up to the past 48 hours).
- Cumulative QC Report – provides information on quality control testing of all panel types. This report does not appear when BD EpiCenter communication is enabled.
- Daily Instrument Report – lists the status of the instrument at the time the report is generated, along with a maintenance checklist.
- Interpretation Rule Set Report – lists the antimicrobial breakpoints of the Interpretation Rule Set defined as the default Rule Set in the Instrument Configuration screen.
- Needs Attention List Report – lists all panels in the instrument's database that have an unignored Needs Attention flag set.
- Organism ID Code List Report – displays all Organism Names, BD and LIS codes for all Organism Names that exist in the BD Phoenix database.
- Panel Lot Report – lists all the panel records for any panel lot number in the current instrument.
- Panel Lot Database Report – lists all the defined panel lots in the current instrument, and provides statistical and reference information on those lots.
- QC Lab Report – lists all Test Strain Organisms that have completed testing and all biochemical and/or antimicrobial MIC results (for a specified Panel Lot #) that exist in the BD Phoenix databases when the report is requested.
- Resident Panel Report – lists the panels contained in stations 1–25 for each tier.
- Lab Report – contains all information for a panel sequence number that exists in the BD Phoenix database, including all information in the Panel Results screen, any special messages, BDxpert Rules that triggered, or Needs Attention Reasons if they exist. This report can only be printed from the Results or Finalization tabs.
- Finalization Summary Report – contains all the panels eligible for finalization at the time the report was requested, as well as finalization status. This report can only be printed from the Finalization tab.
- Panel Inventory Lab Report – prints a Lab Report for all panels in the instrument. This report can only be printed from the Inventory tab.

- Custom Breakpoint Difference Report – prints a listing of differences between old breakpoints and new ones after a Phoenix Update Data or Install/Upgrade operation. This report can only be printed from the Custom Interpretation Rule Set (Configuration) tab.
- Current QC Panel Lot Report – contains information on the most recent QC test for each of the required strains for a panel lot, up to a maximum of 20 strains. This report can only be printed from the Panel Lot Definition tab.
- Historical QC Panel Lot Report – contains information on all tests for a strain for the current instrument, up to 200 tests. This report can only be printed from the Panel Lot Definition tab.

3.8.1 How to Print Reports

To print a report:

- 1 Select the **Reports** tab.
- 2 Highlight the desired report.
- 3 Complete any additional fields (such as an Accession Number for the Accession Lab Report) and select **Print Reports**. (See Section 3.8 Reports Tab for additional requirements.)

Several reports can also be printed from the screens that relate to them (e.g., Needs Attention List Report).

Each of the reports is discussed in greater detail in the sections that follow.

3.8.2 Completed Lab Report

This report contains information for **all panels** whose status became **Complete** during the time selected period (up to the past 48 hours). The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), and Date and Time Printed, Software version / PUD version.

Body of Report: Top Region: Accession #, Isolate #, Sequence #, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Completed Lab Report is sorted by Accession # then by Isolate # within accessions.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any Instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDXpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDXpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report. Any SIR values and Rule # are not reported for QC panels.

ID and AST sections of combination panels print on separate pages of the report.

3.8.3 Accession Lab Report

This report is basically a collection of lab reports for a specified accession number. It provides information for a specified accession number. Such information includes, as applicable: organism ID results, including specific biochemical reactions; AST results including SIR interpretation and MIC; QC pass/fail results; and any BDXpert Rules that were triggered. The report provides the following information:

Header: Report Title, Preliminary indication (if Status is Ongoing, Pending, or Partial Complete; and/or if there are unignored Needs Attention conditions or a Needs Attention condition of Cannot Identify Barcode; and/or if the panel is not eligible for finalization), Laboratory Information (if configured), Date and Time Printed, Software version / PUD version.

Body of Report: Top Region: Accession #, Isolate #, Sequence #, Panel Type, Status (Ongoing, Complete), CRITICAL panel indication, Test Start with time, Test End with time, Instr #/Station (location), Finalized status, Panel Lot # (if QC Lot Support is enabled), and Inoculum Density. The Completed Lab Report is sorted by Accession # then by Isolate # within accessions.

Below this information the organism Final ID is listed. An asterisk next to the Final ID indicates that the ID was changed by the user. Below this, the Media Type (for Yeast ID panels only) appears.

Next any Instrument ID results are listed, along with the Confidence Value for the result. In the lower region of the report, for ID tests, the Biochemical, Instr(ument) Result, and Expected Result are provided. For AST tests, the Antimicrobial, Instr(ument) MIC, Instr(ument) SIR, BDXpert SIR, Final SIR, and Rule # are printed. If any panels have Resistance Markers, BDXpert Rules, Needs Attention flags, or Special Messages, these are printed at the bottom of the report. Any SIR values and Rule # are not reported for QC panels.

ID and AST sections of combination panels print on separate pages of the report.

For QC panels, in addition to the information listed above, the following information is included: Panel Lot # and Expiration Date, Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, and Test Strain. The QC Status of PASS, REVIEW, ERROR, or FAIL is indicated.

3.8.4 Needs Attention List Report

This report lists all the panels in the instrument's database that have an unignored Needs Attention. This report can also be printed from the Needs Attention screen. The report provides the following information:

Header: Report Title, Filtered report notification (if selected), Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Needs Attention Reason, Sequence #, Accession #, Isolate #, Instr #/Station (location), and Status (Ongoing, Complete). If the report is filtered (default selection), an asterisk appears to the left of the Reason for panels with multiple Needs Attention conditions. The report is sorted by the priority of the Needs Attention Reasons (Figure 3-7 – Needs Attention List Report), and by Accession within each Reason type.

Display Panel to Review	Sample ID (Accession)	Accession Number	Isolate Number	Isolate Set
Test Aborted	427290613434		623	
Test Aborted	427290613433		624	
Test Aborted	427290613432		625	
Test Aborted	427290613430		626	

Figure 3-7 – Needs Attention List Report

3.8.5 Resident Panel Report

This report lists the panels contained in stations 1–25 for each tier detected during the last inventory scan. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument #.

Body of Report: Accession #, Isolate #, Sequence #, QC (if panel is QC), Test Start with time, Inoculum Density, Panel Type, Status (Ongoing, Complete), and the highest priority Needs Attention Reason if one exists. The report is sorted by Accession # and then by Isolate # within each Accession.

3.8.6 QC Lab Report

This report lists all QC panels from the Test Start date field entry to the current date. It lists all Test Strain Organisms that have completed testing and all biochemical and/or antimicrobial MIC results (for a specified panel lot number) that exist in the BD Phoenix database. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed, Software version / PUD version.

Body of Report: Panel Lot # and Expiration Date, Test Start and time, Test End and time, Panel Type, Instr #/Station (location), Status (Ongoing, Complete), Tech ID, ID Broth Lot # and Expiration Date, AST Broth Lot # and Expiration Date, Indicator Lot # and Expiration Date, Sequence #, Accession #, Isolate #, Test Strain, Inoculum Density, Media Type (Yeast ID panels only), Instrument ID(s), and QC Status of PASS/FAIL/REPEAT/REVIEW. At the bottom of the report, any Needs Attention reasons or Special Messages are printed. Each Biochemical, along with Instr(ument) Result, and Expected Result are provided, as well as Antimicrobials, Instr(ument) MICs, Expected MICs, and Pass/Fail status.

This report is only available when BD EpiCenter is disabled.

3.8.7 Cumulative QC Report

This report provides information on completed quality control testing of all panel types. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument #.

Body of Report: Selection Criteria: Panel Lot #, Panel Type, and Test Strains selected.

Below this, the panel Sequence # (sort order), QC Status (PASS, FAIL, REVIEW, ERROR), Test Strain, Test Start and Time, Panel Lot #, ID Broth Lot #, AST Broth Lot #, Indicator Lot #, and Tech ID are listed for each panel.

This report is only available when BD EpiCenter is disabled.

3.8.8 Daily Instrument Report

This report lists the status of the instrument at the time the report is generated, and provides areas to record maintenance activities. The Daily Instrument Report can be set to print automatically at a specified time.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, and Software version/PUD version.

Body of Report: Instrument #, Serial #, Instrument Temperature Pass/Fail status, Carousel Rotational Test Pass/Fail status, Power Supply Check Pass/Fail status, Normalizer Panels Sequence #, Pass/Fail status, and Expiration Status (date if expiration is more than 60 days; "expires on date" if expiration is between 60 and 0 days; and "EXPIRED" if the panel is expired) for each tier, and blanks to record the reading, Pass/Fail status, and Tech ID for each of the following maintenance checks: Daily: Instrument Temperature (Status screen), Standard Panel Temperature, Printer Paper Supply: Weekly: Internal Green LEDs, Internal Red LEDs, Internal Amber LEDs, Alert Indicator, and Instrument Audible Alarm. An area is provided for comments at the bottom.

The instrument temperature is considered to have passed when there are no outstanding E01 temperature alerts.

3.8.9 Interpretation Rule Set Report

This report lists the antimicrobial breakpoints of the currently selected Interpretation Rule Set (defined as the default Rule Set in the Instrument Configuration screen). The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Rule Set, Rule Version, and columns for Antimicrobial (sort order), Test Group, Organism Group, Organism Name, S(usceptible) value, and R(esistance) value. Each antimicrobial breakpoint is listed in a separate row of the report.

The Interpretation Rule Set is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

To print reports of ALL the rule sets in the instrument, first select a rule set in System Configuration (Section 3.10.2 System Configuration Sub-Tab), then access **Reports** and select Interpretation Rule Set Report. The currently selected rule set prints. When printing is complete, return to **System Configuration**, and select the next rule set. Then access **Reports** again and print the current Interpretation Rule Set. Continue selecting rule sets and printing until all rule set selections are printed. Remember to return to System Configuration and select the desired rule set to use for interpretations when all printing is complete.

NOTE

Modification of the Interpretation Rule Set should not be performed while there are ongoing panels. This could lead to inaccurate interpretations.

Interpretation Codes:

CLSI or EUCAST

Interpretation Code	Interpretation Name
ACIN_IC	<i>Acinetobacter</i> spp.
AERM_IC	<i>Aeromonas</i> spp.
BURCEP_IC	<i>Burkholderia cepacia</i> complex
ENTC_IC	<i>Enterococcus</i> spp.
ENTERIC_IC	<i>Enterobacteriaceae</i>
NFGNROTH_IC	Nonfermentative GNR, other than ACIN_IC, BURCEP_IC, PSEAER_IC, STEMAL_IC, ACTBACT_IC, CARHOM_IC, EIKCOR_IC
PSEAER_IC	<i>Pseudomonas aeruginosa</i>
STAAUE_IC	<i>Staphylococcus aureus</i>
STAOOTH_IC	<i>Staphylococcus</i> spp., other than STAAUE_IC
STEMAL_IC	<i>Stenotrophomonas maltophilia</i>
STRBET_IC	<i>Streptococcus</i> beta-hemolytic
STROTH_IC	<i>Streptococcus</i> spp., other than STRBET_IC, STRPNE_IC, STRVIR_IC
STRPNE_IC	<i>Streptococcus pneumoniae</i>
STRVIR_IC	<i>Streptococcus viridans</i> group

3.8.10 BDxpert Rule Set Database Report

This report lists each BDxpert rule number and the text describing the rule, whether each rule is enabled/disabled and whether each rule shall trigger automatically/manually in the system. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Rule Set (CLSI, SFM, EUCAST, or Custom), and Based On (CLSI, SFM, or EUCAST if Rule Set is Custom).

Body of Report: Rule #, text of the rule, Enabled/Disabled status, and Automatic/Manual status.

The BDxpert Rule Set Database is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

3.8.11 Organism ID Code List Report

This report prints all Organism Names and Abbreviations for all Organism Names that exist in the BD Phoenix database. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Organism name (sort order), BD Code (abbreviation), and LIS Code (if enabled).

The Organism ID Code List is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

3.8.12 Antimicrobial Code Report

This report prints all antimicrobials and abbreviations for all antimicrobials that exist in the BD Phoenix database from all panel configurations. The report provides the following information:

Header: Report Title, Laboratory Information (if configured), and Date and Time Printed.

Body of Report: Antimicrobial name (sort order), BD Code (abbreviation), and LIS Code (if enabled).

The Antimicrobial Code is a large report. Spooling and printing the report can consume system resources such that other reports cannot be printed until the current one completes.

3.8.13 Lab Report/QC Lab Report

This report (Section 3.4.9 Needs Attention Sub-Tab) contains all information for a panel sequence number that exists in the BD Phoenix database, including all information in the Panel Results screen, any special messages, BDxpert Rules that triggered, or Needs Attention Reasons if they exist. The Lab Report is not accessible from Reports. It can only be printed from Results or Finalization. A Panel Inventory Lab Report prints the same information for all panels listed in the Panel Inventory screen (i.e., resident in the instrument and with final or partial results). Refer to Section 3.4.9 Needs Attention Sub-Tab for additional information.

The QC Lab Report (Section 3.8.6 QC Lab Report) is also accessible from Results. It provides similar information to the Lab Report, but for QC panels.

3.8.14 Finalization Summary Report

This report contains a listing of all the panels eligible for finalization at the time the report was requested, as well as finalization status. The Finalization Summary Report is not accessible from Reports; it can only be printed from Finalization. Refer to Section 3.8.14 Finalization Summary Report for additional information.

3.8.15 Custom Breakpoint Difference Report

This report contains a listing of differences between old breakpoints and new ones after a BD Phoenix Update Data or install/upgrade operation. The Custom Breakpoint Difference Report is not accessible from Reports. It can only be printed from the Custom Interpretation Rule Set (Configuration) tab by selecting the Difference Report button located at the bottom of the screen. Breakpoints that have been customized will not be overwritten with updates from the PUD. The Custom Breakpoint Difference Report will provide the appropriate data to determine if customized breakpoints need to be manually updated to reflect the currently installed PUD. Refer to Table 3-3 – Rules Updates for details on how breakpoints are updated.

3.8.16 Current QC Panel Lot Report

This report contains information on the most recent QC test for each of the required strains for a panel lot, up to a maximum of 20 strains. The report includes information for any instruments whose data has been restored to the current instrument. The Current QC Panel Lot Report cannot be printed from Reports. It can only be printed from Panel Lot Definition. Refer to Section 3.10.9 Panel Lot Definition Sub-Tab.

3.8.17 Historical QC Panel Lot Report

This report contains information on all tests for a strain for the current instrument (only), up to 200 tests. The Historical QC Panel Lot Report cannot be printed from Reports. It can only be printed from Panel Lot Definition. Refer to Section 3.10.9 Panel Lot Definition Sub-Tab.

3.8.18 Panel Lot Report

This report lists all the panel records for any panel lot number in the current instrument. The report first lists clinical panels, then QC panels. Within each of those groups, the report is sorted by Accession # then Isolate #.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Instrument where printed.

Body of Report: Panel Lot # and Panel Type; Accession #, Isolate #, Sequence #, QC (if panel is QC), Test Date, and Status (Pending, Ongoing, Complete) for each panel tested that belongs to the lot.

This report is only available when BD EpiCenter is disabled and QC Lot Support is enabled.

3.8.19 Panel Lot Database Report

This report lists all the defined panel lots in the current instrument, and provides statistical and reference information on those lots.

The report provides the following information:

Header: Report Title, Laboratory Information (if configured), Date and Time Printed, Instrument where printed.

Body of Report: Panel Lot # (sort order, descending); Panel Type; Expiration Date; Extension Date (if Expiration date was extended); Start and End Sequence #s (Range); Definition Date; First and Last Date Used; and number of Panels Used.

This report is only available when BD EpiCenter is disabled and QC Lot Support is enabled.

3.9 Maintenance Tab

Maintenance (see Figure 3-8 – Maintenance Screen) provides several tasks for performing instrument maintenance. There are user tasks for weekly and as needed maintenance. There are other tasks for BD use only.

Under Task Category, select the category of maintenance to be performed:

- Hardware
- Software
- LIS
- Panel
- Field Service
- Internal BD Use


Figure 3-8 – Maintenance Screen

3.9.1 Maintenance Hardware Functions

For detailed information, refer to Section 5.2 Routine Maintenance.

- Test system indicators
- Test Alarm

3.9.2 Maintenance Software Functions

Each of the software functions are described in detail below. Unless otherwise specified, perform the following steps to save data after accessing each software function via the Maintenance tab.

- 1 Insert the USB key, BD part number 443866, into the AIO PC.
- 2 Select the Maintenance screen tab.
- 3 On the Maintenance screen select Software under the Task Category.
- 4 From the Task list select the function to be performed.
- 5 Select Execute; the **Are You Sure?** message is displayed.
- 6 Select OK.

Save System Data Task

Under certain circumstances, BD may advise that system data be saved to a USB key. These circumstances include some error conditions and system malfunctions. The Save System Data function is NOT a backup and cannot be restored by the user.

Save Event Log Task

Save Event Log to Network should only be used when advised by your local BD representative. The Save Event Log to Network option appears only when the instrument is connected to a BD EpiCenter system. This option enables the event log to be saved to an BD EpiCenter system.

Under certain circumstances, a BD representative may advise that the Event Log be saved to the BD EpiCenter system or to a USB key. These circumstances include some error conditions and system malfunctions. The Save Event Log function copies the system event log, which contains logged system messages about various system, instrument, and communications events. The instrument door must be closed, and the instrument must be idle to save the Event Log.

Upgrade Task

Upgrade enables instrument and AIO PC updates to be performed. To perform this operation, the instrument door must be closed. The software update will be recorded on the Daily Instrument Report for reference. To upgrade, follow the step by step instructions in Section 3.9.1 Maintenance Hardware Functions. After OK is Selected, the following window appears:


Figure 3-9 – Instrument Reboot Message

After the upgrade, files on the USB key are validated. The system will copy the upgrade files and reboot.

After the AIO PC finishes rebooting, the following message will be displayed:


Figure 3-10 – Installation Message

The message above will show the status of each step of the upgrade process. The user will not be able to cancel this process once it has started.

When the upgrade completes successfully, the following message is displayed:


Figure 3-11 – Setup Wizard Message

When **Finish** is selected, the system will reboot and the newly installed version of the application on the AIO PC will synchronize with the instrument. If the upgrade includes new software for internal components of the instrument, these updates will take place immediately and the use of the application will be temporarily blocked. Once the instrument software has been successfully updated, the application will be available for use.

If custom breakpoints are in place, print the Custom Breakpoint Difference Report (See 3.10.6 Custom Interpretation Rule Set Sub-Tab)

NOTE

When the software is updated, or when the PUD install is performed, if any of the standard interpretation rule sets have changed, new rule sets are installed into the system database. Since custom rule sets are based on standard rule sets, a custom interpretation rule set will be merged with new standard rules. If custom breakpoints are used, be sure to print out the Custom Breakpoint Difference Report (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab) after each system update or PUD install.

Save User Data Task

Save User Data allows the back up of Configuration parameters to a USB key. This includes Custom Interpretation Rule Set Configuration and BDxpert Rules Configuration. The information that is saved is for the current instrument only. However, information saved at one instrument can be restored on another instrument. It is recommended that Configuration parameters be stored on a USB key in the event of a system failure.

Restore User Data Task

Restore User Data enables the restoration of the saved configuration parameters. This data includes: Custom Interpretation Rule Set Configuration, and BDxpert Rules Configuration. To restore, the instrument door must be closed.

BD Phoenix Update Data (PUD)

A BD Phoenix Update Data may be provided from time to time to update BD Phoenix M50 instrument databases and support files. These updates do not affect or change the basic instrument application software. For this reason, the Status screen shows both the software version and the PUD version near the top of the screen.

BD Phoenix Update Data (PUD) enables the update of numerous data files in the instrument, such as antimicrobial breakpoints and rules, QC data, drugs, organisms, etc.

NOTE

After a PUD upgrade, a Custom Breakpoint Difference Report should be printed and reviewed to determine if manual updates are required for custom breakpoints. Refer to 3.8.15 Custom Breakpoint Difference Report for additional information on this report. Refer to 3.10.6 Custom Interpretation Rule Set Sub-Tab for additional information on how custom interpretation rules are updated with a PUD upgrade.

Backup SQL Database Task

This feature provides data to BD Service for instrument troubleshooting purposes.

3.9.3 Maintenance LIS Functions

Each of these functions are described in detail below. Unless otherwise specified, perform the following steps to save data after accessing each software function via the Maintenance tab.

- 1 Insert the USB key, BD catalog number 443866 into the AIO PC.
- 2 Select the Maintenance screen tab.
- 3 On the Maintenance screen select LIS under the Task Category.
- 4 From the Task list select the function to be performed.
- 5 Select Execute; the **Are You Sure?** message is displayed.
- 6 Select OK.

Save LIS Codes Task

Save LIS Codes enables all the Organism and Antimicrobial LIS codes you have defined to be saved. This produces a text file that can be edited on a PC, which might be quicker for some users who have many edits to perform. Then, the edited codes can be restored back to the BD Phoenix M50 instrument. The function also enables codes to be copied from one instrument to another. Save LIS Codes appears on the Maintenance tab only if LIS Communications is enabled.

Entries in the text field consist of:

Identifier (ORG, QC_ORG, DRUG) | BD Code | Short Name/Drug | LIS Code

Codes cannot be added or deleted, and only the LIS Code portion can be modified. If another field is changed, the instrument will not restore the codes.

Restore LIS Codes Task

Restore LIS Codes enables the restoration of the Organism and/or Antimicrobial codes that were previously saved. The restore operation completely overwrites the existing Organism/Antimicrobial LIS Code database. Codes will not be restored if any field other than LIS Code was modified, or if LIS Code was entered in an incorrect format. If this happens, an error log is written, and the file can be reviewed to see what caused the error. Restore LIS Codes appears on the Maintenance menu only if LIS Communications is enabled.

3.9.4 Maintenance Panel Functions

Each of these functions are described in detail below. Unless otherwise specified, perform the following steps to save data after accessing each software function via the Maintenance tab.

- 1 Insert the USB key, BD catalog number 443866 into the AIO PC.
- 2 Select the Maintenance screen tab.
- 3 On the Maintenance screen select Panel under the Task Category.
- 4 From the Task list select the function to be performed.
- 5 Select Execute; the **Are You Sure?** message is displayed.
- 6 Select OK.

Save Panel Lot Definitions Task

This enables the transfer of lot definitions (and QC panel results) to other instruments so that the records can be viewed/used there. Save Panel Lot Definitions saves the defined panel and related data (Sequence #s, Expiration Dates, etc.), and QC results for any strains tested in those. Save Panel Lot Definitions only applies when QC Lot Support is enabled.

Restore Panel Lot Definitions Task

This makes Panel Lot definitions and QC panel results transferable in labs that use multiple instruments, so that a lot only has to be defined once using the box (carton) label. Restore Panel Lot Definitions only applies when QC Lot Support is enabled.

Save Panel Configuration Task

Save Panel Configuration enables the backup of the instrument's panel configuration. The information saved is for the current instrument only. However, information saved at one instrument can be installed at another instrument. This enables the presence of consistent panel configurations among all the instruments.

Install Panel Configuration Task

Install Panel Configuration enables the update of the instrument's panel configurations. To perform this operation, the instrument door must be closed.

NOTE

This operation should not be performed while there are ongoing panels.

3.9.5 Maintenance Field Service Functions

For Field Service use only.

3.9.6 Maintenance Internal BD Use Functions

For internal BD use only.

3.9.7 Maintenance Event Log Tab

The Event Log provides a list of messages generated by LIS communications. Note that the screen (see Figure 3-12 – Event Log Screen) will not contain messages if the BD EpiCenter system is connected. The messages represent status messages that have occurred during LIS communications such as: query messages, log entries, and interface messages.

The Event Log list of LIS messages may be filtered based on date range. There is also a Find function to search for specific message content.


Figure 3-12 – Event Log Screen

3.10 Configuration Tab

To access the configuration functions, select the **Configuration** tab. The configuration screen then appears with sub-tabs for each specific configuration: Users, System, Communications, Custom Interpretation Rule Set, BDXpert Rules, Rapid Reporting, and Instrument. When the desired configuration parameters have been entered or modified, select **Save** to make the changes permanent.

* Note that the Custom Interpretation Rule Set and BDXpert Rules Configuration tabs do not appear if the BD Phoenix M50 instrument is connected to and communicating with an BD EpiCenter system.

3.10.1 Users Sub-Tab

Refer to Section 2.4 User Access and Management Summary.

3.10.2 System Configuration Sub-Tab

The following parameters can be set in System Configuration (see Figure 3-13 – Systems Configuration Screen):

GENERAL System Settings

System Number

Select the system identification number. The default setting is 0. Choose a number from 1 to 99. If there is only one instrument, leave this value set at 1.

Rule Set

Select the rule set that the interpretation engine is to use. Only one rule set can be used. Choose from the following selections:

CLSI (Clinical and Laboratory Standards Institute)

EUCAST (European Committee on Antimicrobial Susceptibility Testing)

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Custom (defined in Section 3.10.6 Custom Interpretation Rule Set Sub-Tab)

Rule Version

This read-only field shows the current version of the rule set selected in the previous field.

Alarm (audible) Volume

Select the volume of the instrument's audible alarm. The default setting is 5. Select from 0 (audible alarm off) to 10 (loudest). Only the volume of Alert and Activity alarms (see Section 2.3.2.8 Audible Tones and Alarms) is affected by this setting.

Resistance Marker Notification

This field appears only when the BD Phoenix M50 instrument is connected to a BD EpiCenter system but is not communicating with it. Enable the instrument to operate in a standalone mode when communications with BD EpiCenter is interrupted. When this field is enabled, the BDxpert System Active field appears.

BDxpert System Active

This field activates the entire BDxpert System rules which includes CLSI, EUCAST, SFM, or Custom Interpretation Rules.

This field does not appear if the BD EpiCenter system is attached and communicating. However, if communication with the BD EpiCenter system is lost, the field reappears. This enables the activation of BDxpert rules interpretations in the standalone BD Phoenix M50 instrument. When communications with the BD EpiCenter system is restored, BDxpert rules interpretations are once again performed at the BD EpiCenter system, and this field is removed from the screen.

NOTE

Disabling all BDxpert rules also disables detection of Resistance Markers (e.g., ESBL), except those triggered by 1500-series rules.

Rapid Completion

This field enables the instrument to provide BDxpert AST results (SIR) before determining actual MIC values. The instrument MIC values are provided as soon as they can be accurately determined. Within a test panel, some MIC values may be available earlier than others. The rapid completion feature can be used to predict resistance for uncompleted antibiotics using the ID alone (intrinsic resistance), or ID with completed MICs of related antibiotics, or resistance marker tests (BL, ESBL). The BDxpert system is used to make these predictions. This can be useful in situations where, for example, the results for drugs that have not yet received MICs would be of no clinical value based on the other results that are already available. Antimicrobials with Rapid Complete BDxpert interpretations are indicated by a **C** in the MIC column on Results screens and Lab reports.

If Rapid Completion is selected at the BD Phoenix M50 instrument and the instrument is connected to BD EpiCenter, then it must also be enabled (checked) at the BD EpiCenter.

OPTIONS

Abbreviated Lab Report

Enabling this field causes the system to print a shortened version of the Lab Report. The shortened version does not contain results for the biochemical (ID) reactions. Disable this field to print the standard full-length Lab Report. The default is disabled (full length). Note that QC Reports always print standard full-length Lab Reports.

QC Lot Support

This field enables the QC Lot feature, which can be used to facilitate panel lot QC testing and tracking (see Section 3.4.6 Lot Information Sub-Tab).

Daily Instrument Report Printing

System Check Box

Enabling this field causes the Daily Print Hour and Daily Print Minute fields to appear, and tells the instrument to print a Daily Instrument Report automatically at the time specified in these fields. If automatic printing is disabled, the report can still be printed at any time from Reports.)

Laboratory Information

Information entered in the Laboratory Information window prints in the header of system reports.

Name

Enter the laboratory name, up to 40 characters.

Address

Enter the address, up to 40 characters.

City

Enter the city (as well as state and zip code, if desired), up to 40 characters.

Director

Enter the name of the laboratory director, up to 40 characters.

Inoculum Density

Inoculum Density enables the default McFarland concentration for inoculum to be set. Select 0.25 (acceptable density of 0.20–0.30) or 0.5 (acceptable density of 0.50–0.60) for Gram Positive and Gram Negative panels (Strep panels use only 0.5). The default density is 0.5 for Gram Negative and Gram Positive panel types. The density for Yeast ID panels is fixed at 2.0 McFarland (acceptable density of 2.00–2.40). Inoculum density is applicable only to panels that have an ID side.

Default Media Type

Yeast – This field enables the selection of a default media type that appears during Panel Login when a Yeast ID panel sequence number is scanned or typed in. From the drop down box, highlight the desired media and select save to store the configuration setting. A default media type does not have to be specified, but a media type must be selected when logging in Yeast ID panels.


Figure 3-13 – Systems Configuration Screen

3.10.3 Communications Sub-Tab


Figure 3-14 – Communications Configuration Screen

The Communications Configuration screen (Figure 3-14 – Communications Configuration Screen) enables/disables and configures communications for the BD Phoenix M50 instrument with a compatible LIS (Laboratory Information System). It also enables BD representatives to enable/disable/adjust communications with the BD EpiCenter advanced data management system. Only LIS or BD EpiCenter communications can be enabled. If there is a BD EpiCenter system connected and communications with a LIS system is required, BD EpiCenter can be configured to communicate with the LIS.

LIS operations are discussed in Section 2.8 LIS Operations. The AIO PC displays the LIS connection status icon the upper left side of status screen.

LIS Enabled

Touch the **LIS Enabled** field to establish a connection. When LIS Enabled is checked, the following fields appear: Network Configuration, Options, and Results Upload Options.

Network Configuration

Baud

Available choices are: 2400, 4800, 9600 (default), 14400, 19200, 38400.

Data Bits

Available choices are: 7, 8 (default).

Parity

Available choices are: None (default), odd, even.

Stop Bits

Available choices are: 1 (default), 2.

Packed Frames

Select whether packed frames can be used for serial communications with the LIS system. Enable this field to allow the BD Phoenix M50 instrument to send multiple records per frame. A disabled field indicates that one record per frame is uploaded to the LIS.

Options

Send Interpretation Results

Enabling this option causes the final SIR values for antimicrobials to be included in the Results record uploaded to the LIS. Its default value is enabled.

Unsolicited Queries

Enabling this option causes the BD Phoenix M50 instrument to request panel information from the LIS if the panel is placed into the instrument lacking an organism ID for AST panels or Combination panels with only the AST portion of the panel enabled.

This field's default value is disabled, which means the BD Phoenix M50 instrument will NOT request missing information from the LIS.

Send When Placed in Instrument

Enabling this option causes the BD Phoenix M50 instrument to send a Results upload to the LIS when the panel is placed into the instrument. This field's default value is enabled.

ASTM Byte Mode Comments

This field only appears when Japanese is selected for the language and causes the ASTM standard to be followed (when enabled) or not followed (when disabled). When the ASTM standard is not followed, escape sequence characters that are usually rejected are instead accepted.

This field's default value is enabled.

Results Upload Options

Only one option from the list below can be selected.

Solicited

Enabling this option causes the BD Phoenix M50 instrument to upload Results records only when the LIS requests the information. Its default value is disabled.

QC panels and orphan panels are uploaded only when solicited by the LIS.

Send on Finalization

Enabling this option causes the BD Phoenix M50 instrument to upload Results records only when the panel is finalized. Its default value is disabled.

Send on Completion

Enabling this option causes the BD Phoenix M50 instrument to upload Results records only when the panel status becomes complete or when a change is made to a complete panel. Its default value is enabled.

Send as Available

Enabling this option causes the BD Phoenix M50 instrument to upload Results records at the following times: when an ID is determined; when an AST result is determined; when there is a change to a panel record that already has at least partial results. Its default value is disabled.

Send at Fixed Time

Enabling this option causes the BD Phoenix M50 instrument to upload Results records for a panel with partial results at a fixed, specified time. The default value is disabled.

NOTE

Send as Available and Send at Fixed Time results should be considered PRELIMINARY until the panel is complete.

Results records are uploaded from the BD Phoenix M50 instrument to the LIS system. These records are covered in detail in Section 2.8.2 Results Upload Records.

3.10.4 Organism Configuration Sub-Tab

The Organism Configuration screen allows the LIS codes for the organisms in the BD Phoenix M50 instrument database to be edited. LIS codes for organisms must be unique.

The Organism Configuration screen enables the selection of the organism (or Test Strain for QC organisms) to be edited.

To edit an organism LIS code:

In the Organism Configuration screen, select the name in the list and type in the new LIS Code in the Modify Organism ID LIS Code area at the bottom of the screen. Select Save LIS Code to save the changes. When all modifications to the LIS codes have been completed, save and exit the Communications Configuration screen.

3.10.5 Antimicrobial Configuration Sub-Tab

The Antimicrobial Configuration screen allows the LIS codes for the antimicrobials in the BD Phoenix M50 instrument database to be edited. LIS codes for antimicrobials must be unique.

This screen also enables the selection of the antimicrobial to be edited. Antimicrobials are presented alphabetically by antimicrobial name.

To edit an antimicrobial LIS code:

In the Antimicrobial Configuration screen, select the antimicrobial and type in the new LIS code in the Modify LIS Code window at the bottom of the list. Select Save LIS Code to save the modifications. When all modifications to the LIS codes are complete, save and exit the Communications Configuration screen.

3.10.6 Custom Interpretation Rule Set Sub-Tab

Custom Interpretation Rule Set Configuration enables interpretation rules to be tailored to the specific needs or requirements of your laboratory. The current rule set (whether default or already customized) is used as the basis for custom rules. The current interpretation rule set from Reports can be printed (see 3.8 Reports Tab). Note that this configuration function is not available if the BD Phoenix M50 instrument is connected to and communicating with a BD EpiCenter system.

Periodically, a software update or BD Phoenix Update Data (PUD) operation may install new breakpoints. When a software or PUD update has been performed, print a Custom Breakpoint Difference Report by selecting the Difference Report button at the bottom of the screen. Breakpoints that have been customized will not be overwritten with updates from the PUD. The Custom Breakpoint Difference Report will provide the appropriate data to determine if customized breakpoints need to be manually updated to reflect the currently installed PUD. Refer to Table 3-1 – MIC/SIR Values and Causes for details on how breakpoints are updated. Refer to Section 3.8.15 Custom Breakpoint Difference Report for additional details on this report.


Figure 3-15 – Custom Interpretation Rule Set Screen

Rules (Standards) Updates

The chart below depicts the logic of updates to Rules Sets. The top of the chart shows the effect of rules changes when a customized rule set is present. The bottom of the chart shows the rest of the cases for rules updates.

In the chart, the letters A, B, and C represent the entirety of a rule set (standard). For example, CLSI Rule Version M100_S25 might be indicated by value A. An update, delivered via PUD, might install CLSI Rule Version M100_S26, which might be indicated by the value B. The customization of A might be represented by C. Therefore, by following the chart below, the effect of the update is that value C (the customized rules based on Version M100_S25) would be retained.

	Old Standard	New Standard	Old Custom	Result after PUD (New Custom)
Custom rules are present	A	A	B	B
	A	B	C	C
	–	–	A	A
	–	A	A	A
	–	A	B	A
	A	–	B	B
No custom rules	A	A	A	A
	A	B	A	B
	–	A	–	A
	A	A	–	A
	A	B	–	B
	A	–	A	–

Table 3-3 – Rules Updates

NOTE

An interpretation rule is considered custom if any of the following are modified: Susceptible Breakpoint, Resistant Breakpoint, DTG.

To use the custom rule set for interpretation, go to the **System Configuration** tab and select **custom** from the Rule Set drop down menu.

Custom Interpretation Rule Set – Based On:

The Based On drop down list selects the standard rule upon which the custom rules are to be based. The standard rule sets are: CLSI, EUCAST, and SFM.

NOTE

If a rule set other than the current rule set is selected, the currently defined custom interpretation rules are deleted.

Custom Rule Set Table

The main window on this screen shows the current rules in the rule set. The following values are shown from left to right (items are explained below under To modify a rule):

Check box next to a rule (indicates the rule is marked for deletion)

Antimicrobial

DTG (Drug Test Group - see information below)

Org (anism) Group

Organism

S (usceptible Value)

R (esistant Value)

To add a rule:

- 1 Select **Add**. A new rule is added to the end of the list and selected for modification.
- 2 Modify the rule attributes in the Selected Antimicrobial Name Details window and then save.
- 3 Once all new rules have been added select **Save All** to save the new rules and exit the custom interpretation rule set configuration screen.

To modify a rule:

- 1 Select the row containing the rule to be modified. The row will be highlighted and the rule attributes will be populated in Selected Antimicrobial Name Details. Note - Do not check the field next to the rule. This is used to delete rules.
- 2 Modify the rules attributes and Save when all modifications are complete.

DTGs are derived from the AST Standards and are used to categorize antibiotics into distinct groups. For each drug, the groups are specific to the organism group and the recommended utilization for that drug.

DTGs have no significance in the BD Phoenix system alone. They are used in conjunction with the BDxpert rules and are necessary for interface with the BD EpiCenter system.

Generally, these are divided into seven groups shown below. Only the A, B, C and U codes may be reported on the Lab report(s). B and C groups are only reported when they are promoted to A by the BDxpert system.

- A Always tested and always reported.
- B Usually tested, but not always reported.

- C Sometimes tested and not always reported.
- U Urinary tract specific drug. Tested and reported for source urinary tract.
- O Other drugs that may be tested. Will not be reported unless changed by user.
- I Investigational drugs. Drugs not approved for clinical use are never reported.
- N Not grouped by standard. Will not be reported unless changed by user.

These drug testing codes are used for two purposes in the BD Phoenix M50 instrument. First, drug testing codes can be altered to more closely match the antibiotic formulary and drug utilization guidelines within an institution. The initial codes will be determined by the AST interpretive standard that have been selected. Allowable changes are shown below.

Starting Code	Allowable (Recommended) Changes
A	B, C, U, O
B	A, C, U, O
C	A, B, U, O
U	A, B, C, O
O	A, B, C, U
I	None
N	A, B, C, U

The second application of the drug testing codes is for the promotion or suppression of drug results for the chartable report. This application is driven by the rules in the BDxpert system. If the rules alter the drug testing codes, these will be reflected on the Lab Report generated by the BD Phoenix M50 instrument. The promotion and suppression actions are shown below. At the BD EpiCenter level, the drug testing codes will be used to determine which drug results actually appear on the chartable report.

Starting Code	Promotion	Suppression
A	Not applicable	C
B	A	C
C	A	Not applicable
U	Not applicable	C
O	Not allowed	Not allowed
I	Not allowed	Not allowed

When the desired rules have been modified, select **Save All** to save modifications to the rule set and exit Custom Interpretation Rule Set Configuration.

To delete a rule:

- 1 Check the field next to the rule to be deleted.
- 2 Select **Delete**.
- 3 Select **OK** in the Are you sure message box.
- 4 **Save All** to save the modifications and exit the Custom interpretation rule set screen.

3.10.7 BDxpert Rules Sub-Tab

BDxpert Rule Configuration allows individual BDxpert rules to be enabled or disabled, or triggered automatically or manually according to the specific needs or requirements of each laboratory. The existing BDxpert rule set can be printed by selecting BDxpert Rule Set Database Report from Section 3.8 Reports Tab. Note that this configuration function is not available if the BD Phoenix M50 Instrument is connected to, and communicating with an BD EpiCenter system. See Figure 3-16 – BDxpert Rule Configuration Screen.

BDxpert Rule Set fields:

Rule Set

The rule set currently being displayed is shown in this field. Values are CLSI, SFM, EUCAST, or Custom. Only rules in the selected set appear in the Rule Set field.

Based On

If the rule set is Custom, the rule set on which custom values are based is shown in this field.

Rule Description

The Rule Description field shows the text of the rule selected.

Figure 3-16 – BDxpert Rule Configuration Screen

Rule Set

This screen presents a summary of the status of each of the BDXPert rules and enables rules to be selected for modification.

To enable or disable a rule, Select **Enable/Disable** to toggle the status. When disabled, a rule is not available for use in the system.

To set the rule to automatic or manual trigger, select **Automatic/Manual** to toggle the status. When set to Automatic, a rule will be triggered automatically by the system. When set to manual, the rule must be accepted by the user before it is applied to results.

While 1500-series rules can be selected for viewing, they are fixed as enabled and automatic. These settings cannot be changed. When all parameters have been modified, select **Save**.

3.10.8 Rapid Reporting Sub-Tab

Rapid Reporting Configuration specifies the criteria for notification of critical panels and resistance markers. Notifications for critical panels can originate from an automatic Lab Report printing, an audible alarm, and/or by an upload to the LIS when certain types of results are obtained in non-QC panels. The types of results, which are selected, are ID only, partial results, and complete results. For Resistance Markers, designate the type of notification desired when resistance markers (that are enabled) are triggered. For both, disabling both printing and audible alarms can be chosen. In this case, the only notification is that the panel is shown in red text on the Panel Inventory screen.

Rapid reporting alarms exist in two states: unacknowledged and acknowledged. Unacknowledged alarm panels continue to display a red ! (exclamation mark) on the Inventory screen. Acknowledgment is achieved when an automatic Lab Report for the panel prints or when the panel is viewed in Results. Audible alarms continue to sound for unacknowledged alarms until the **SILENCE ALARM** button is pressed.

Critical Panel Notification

The Critical Panel Notification field allows parameters to be set for notification on critical panels. Critical panels are any panels where special notification occurs when results become available. Notification can be audible and/or by immediate printing of a Lab Report (or neither). Panels can be marked as critical when they are logged in via Panel Login, or subsequently on the Results screen.

Auto Print

When enabled, the system prints a Lab Report whenever the selected results parameters occur (ID Only, Partial, or Complete – see below). This field is disabled by default (auto print off). Both Auto Print and Audible Alarm may be enabled simultaneously.

Upload to LIS

This only appears if LIS Communications is enabled. The default value is disabled. When enabled, the system uploads results records whenever the selected results parameters occur (ID Only, Partial, or Complete – see below). This field overrides the results upload values set in Communications Configuration.

Audible Alarm

When enabled, the system sounds an audible alarm whenever the selected results parameters occur (ID Only, Partial, or Complete – see below). This field is enabled by default (audible alarm on). Both Auto Print AND Audible Alarm may be enabled simultaneously.

Tone Select

This field allows the selection of the tones sounded by the critical panel audible alarm. It appears only if Audible Alarm is checked (enabled). This field is set to Tone 1 by default. The system sounds a sample tone each time the field is selected or its contents are changed.

NOTE

Enable one of the following three fields in order to receive any critical panel notifications.

ID Only

When enabled, critical panel notification is performed only when an organism ID is detected on a panel. This field is set to disabled by default (alarm/report on partial results). If this field is enabled, Partial Results and Complete cannot be enabled.

Partial Results

When enabled, critical panel notification is performed when partial panel results are obtained. This field is set to enabled by default (alarm/report on partial results). If this field is enabled, ID Only and Complete cannot be enabled.

Complete

When enabled, critical panel notification is performed only when complete panel results are obtained. This field is set to disabled by default (alarm/report on partial results). If this field is enabled, ID Only and Partial Results cannot be enabled.

Resistance Marker Notification

This field does not appear if the BD EpiCenter system is connected.

Auto Print

When enabled, the system prints a Lab Report whenever a resistance marker is triggered. This field is set to disabled by default.

Upload to LIS

This appears only if LIS Communications is enabled. The default value is disabled. When enabled, the system uploads results records whenever the selected results parameters occur (ID Only, Partial, or Complete – see above).

This field overrides the results upload values set in Communications Configuration (Section 3.10.3 Communications Sub-Tab).

NOTE

Enable one of the following two fields in order to receive any resistance marker notifications.

All Panels

When enabled, resistance marker alarms are generated for all panels in which markers occur. This field is disabled by default. If this field is enabled, Critical Panels Only cannot be enabled.

Critical Panels Only

When enabled, resistance marker alarms are generated only for critical panels in which markers occur. This field is enabled by default. If this field is enabled, All panels cannot be enabled.

Resistance Marker

This field enables or disables the notification of individual resistance markers. The field shows the abbreviation and text of the Marker, as well as whether it is currently enabled or disabled. By default, all resistance markers are enabled. Resistance marker alarm notification occurs only for enabled markers.

3.10.9 Panel Lot Definition Sub-Tab

Panel Lot Definition enables the definition of panel lots to facilitate QC testing and tracking. Panel lot definition saves the lot number, beginning and ending panel sequence numbers, definition date, and first and last used dates for a lot. Current and historical QC results can be viewed and printed. Note that this configuration function is not available if the BD Phoenix M50 instrument is connected to and communicating with an BD EpiCenter system.

This tab's activities are:

- Define new panel lots
- Recall existing panel lots
- Review current and historical QC results for panel lots
- Print Panel Lot (current, historical) reports

When a panel lot definition is saved, any Test Strain recommendations that exist are displayed in the Current QC Results window. There is a toggle feature between Current and Historical QC Results windows. Current results show a summary of the most recent test for each strain for all instruments, up to 20 entries. Historical results show all tests for a strain for the current instrument only, up to 200 entries. In order to display results from other instruments, Panel Lots must be saved at those instruments and restored to the current instrument (see Restore Panel Lot Definitions Task in Section 3.9.4 Maintenance Panel Functions).

Field contents of these windows are explained below.

Panel Lot Number

This field shows the lot number for the panel. To enter a new panel lot number, scan the lot number barcode (next to the L symbol) on the box. The panel lot number is encoded in the barcode. The lot number cannot be typed in this field when defining new panel lots.

To recall a panel lot definition, select a panel lot from the drop down list, or type the first few numbers of the panel lot number to jump to that portion of the list, or scan the panel lot barcode.

Expiration Date

This field shows the lot's expiration date. This information is encoded in the panel lot barcode and is filled in automatically when the panel lot barcode is scanned. The expiration date cannot be typed in when defining new panel lots. This date can only be changed by scanning a new panel lot barcode provided by the BD Quality department (only if an administrator is logged in).

Start Sequence Number

This field shows the starting barcode sequence number for the panel lot. To enter a new Start Sequence number, scan the lowest sequence number in the lot (next to the symbol on the box). This field can only be completed by scanning the panel sequence barcode; the sequence number cannot be typed in this field.

End Sequence Number

This field shows the ending barcode sequence number for the panel lot. To enter a new End Sequence number, scan the highest sequence number in the lot (next to the symbol on the box). This field can only be completed by scanning the panel sequence barcode; the sequence number cannot be typed in this field.

Definition Date

This field shows the date on which the panel lot was defined. It is completed automatically when the panel lot definition is saved. This is a read-only field that cannot be changed by the user.

First Used Date

This field shows the first date on which a panel from this lot was used. It is completed automatically when the panel is placed in the instrument. This is a read-only field that cannot be changed by the user.

Last Used Date

This field shows the last date on which a panel from this lot was used. It is completed automatically when the panel is placed in the instrument. This is a read-only field that cannot be changed by the user.

Where Defined

This field appears only if the panel lot was defined at a different number instrument than the current one. It shows the instrument number where the panel lot was defined. Panel lot definitions can be saved at one instrument and restored to another instrument through the Save/Restore Panel Lot Definitions functions found on the Maintenance Tab (see Section 3.9 Maintenance Tab).

Current QC Results

This window appears only after a panel lot is saved. When a panel lot definition is saved, any QC Test Strain recommendations that exist are displayed in the Current QC Results window. Only the most recent QC test is shown for each test strain, up to a maximum of 20 tests. Results for other instruments can be shown in this window by saving those results at that instrument and restoring them at the current one (see Save Panel Configuration Task and Restore Panel Lot Definitions Task in Section 3.9 Maintenance Tab).

The following read-only fields are shown in the Current Results window: Test Strain (sort order by name); Test Date; Sequence Number; (Test) Type (ID, AST, or ID/AST); QC Status (None for uncompleted test; Pass, Fail, Review, or Error for completed tests).

A Current QC Panel Lot Report is available by exiting the Panel Lot Definition tab and selecting the Reports tab. In addition to the fields listed above, the report includes standard report header information (Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument number where the report was printed); Panel Lot number; Panel Type; Expiration Date; Instrument number for each strain tested (blank if instrument is the same); and Tech ID.

Historical QC Results

This appears only after a panel lot is saved. Historical results show all tests for a strain for the current instrument only, up to 200 entries. When the desired panel is highlighted, move to the QC Panel Results screen, or print a Historical QC Panel Lot Report.

The following read-only fields are shown in the Historical Results window: Test Strain (drop down selection box); Test Date (sort order); Sequence #; (Test) Type (ID, AST, or ID/AST); QC Status (None, Pass, Fail, Review, Error); Tech ID.

A Historical QC Panel Lot Report is available by exiting the Panel Lot Definition tab and selecting the Reports tab. In addition to the fields listed above, the report includes standard report header information (Report Title, Laboratory Information (if configured), Date and Time Printed, and Instrument number where the report was printed); Panel Lot number; Panel Type; and Expiration Date.

3.10.10 Instrument Sub-Tab

BD field service only.

4 – Routine Panel Operation

4.1 General

This section describes routine panel operations: preparation, quality control and automatic association.

4.2 Storage and Handling

BD Phoenix Panels

Panels are individually pouched and packaged in a box of 25. Panels must be stored unopened at room temperature (15–25 °C). Do not refrigerate or freeze. Visually inspect the packaging for holes or cracks in the foil package. Do not use if the packaging appears to be damaged. Do not use the panel if there is no desiccant or if the desiccant pouch is torn. If stored as recommended, the panels will retain expected reactivity until the date of expiration.

BD Phoenix ID Broth

Tubes are packaged as 100 tube packs. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store BD Phoenix ID Broth tubes at 2–25 °C. Expiration dating is shown on the tube label.

BD Phoenix AST Broth and BD Phoenix AST-S Broth

Tubes are packaged as 100 tube packs. Visually inspect the tubes for cracks, leaks, etc. Do not use if there appears to be a leak, tube or cap damage or visual evidence of contamination (i.e., haziness, turbidity). Store BD Phoenix AST Broth tubes at 2–25 °C. Expiration dating is shown on the tube label.

BD Phoenix AST Indicator Solution and BD Phoenix AST-S Indicator Solution

The indicator solution is individually pouched and packaged as a package of 10 dropper bottles. Visually inspect the bottle for cracks and/or leaks. Do not use if there appears to be a leak, bottle or cap damage, or any change from a dark blue color. Store the BD Phoenix AST Indicator Solution at 2–8 °C. Each bottle contains enough solution to test up to 100 panels. Expiration dating is shown on the box, pouch, and bottle label and is for unopened bottles. An opened bottle will be stable for up to 14 days if stored at 2–8 °C. Be sure the bottle is held vertically when dispensing the AST Indicator Solution.

4.3 Preparing Panels

The BD Phoenix system is not for use directly with clinical specimens. Only pure culture isolates of aerobic and/or facultatively anaerobic gram-negative, gram-positive, and yeast organisms are acceptable for testing. The test isolate must be a pure culture. It is recommended that cultures be 18 to 24 hours old for gram-negative and gram-positive organisms and 18 to 48 hours old for yeast organisms. For AST testing in the BD Phoenix system, isolates recovered from non-selective media are recommended. It is recommended that media containing antibiotics not be used for organisms to be tested in the BD Phoenix system except those specifically claimed in the table below. Selective

media may inhibit some strains of bacteria and yeast, therefore caution must be used when selecting isolated colonies from these media. Use isolates from a blood agar plate such as BD Trypticase Soy Agar with 5% Sheep Blood. The Media selected during login of Yeast ID panels refers to the media on which the organism was grown.

Other recommended media are included in the table below:

Recommended Media	Approved Use			
	ID	AST	Strep	Yeast ID
BD BBL™ CHROMagar™ Orientation	Yes	Yes ¹	X	X
Bromthymol Blue (BTB) Lactose Agar	Yes ⁴	Yes	X	X
Chocolate Agar	Yes	Yes	Yes ²	Yes
Columbia Agar with 5% Horse Blood	Yes	Yes	Yes ³	X
Columbia Agar with 5% Sheep Blood	Yes	Yes	Yes	Yes
Columbia CNA Agar with 5% Sheep Blood (Gram Positives)	Yes	X	Yes	X
Cystine-Lactose-Electrolyte-Deficient (CLED) Agar	Yes ⁵	Yes	X	X
Dey/Engley (D/E) Neutralizing Agar (Gram Negatives)	Yes	X	X	X
Eosin Methylene Blue (Gram Negatives)	Yes	Yes	X	X
Hektoen Enteric Agar (Gram Negatives)	Yes	X	X	X
MacConkey Agar (Gram Negatives)	Yes	Yes	X	X
Phenylethyl Alcohol Agar (Gram Positives)	Yes	X	Yes	X
Sabouraud Brain Heart Infusion Agar - SAB HI (Yeast)	X	X	X	Yes
Sabouraud Dextrose Agar (Yeast)	X	X	X	Yes
Sabouraud Dextrose Agar-Emmons (Yeast)	X	X	X	Yes
BD Trypticase Soy Agar with 5% Sheep Blood	Yes	Yes	Yes	Yes
BD Trypticase Soy Agar with Lecithin and Tween 80	Yes	X	X	X
BD Trypticase Soy Agar without Blood	Yes	X	X	X
Xylose Lysine Desoxycholate Agar (Gram Negatives)	Yes	X	X	X

¹ The use of CHROMagar Orientation may produce false susceptibility results when testing erythromycin with Gram positive organisms. Antimicrobial susceptibility test results should be confirmed using BD Trypticase Soy Agar with 5% Sheep Blood.

² This media type should not be used for Streptococcal identification with SMIC/ID panels. Chocolate Agar may be used for Streptococcal susceptibility testing only.

³ The use of Columbia Agar with 5% Horse Blood may produce significantly higher MIC for SXT with *Streptococcus* species, which may result in false resistance. Antimicrobial susceptibility test results should be confirmed using BD Trypticase Soy Agar with 5% Sheep Blood.

⁴ The use of Bromthymol Blue Lactose Agar with Gram Positive organisms should be restricted to *Staphylococci* for both the 0.5 and 0.25 GP systems.

⁵ The use of Cystine-Lactose-Electrolyte-Deficient Agar with Gram Positive organisms should be restricted to *Staphylococci* for the 0.25 GP system.


Table 4-1 – Recommended Media

It is recommended that the applicator swabs are sterile cotton swabs. Polyester swabs are not recommended. The quality of applicator swabs may vary from vendor and on occasion, loose fibers may dislodge from the swab affecting McFarland readings.

The usefulness of the BD Phoenix system or any other diagnostic procedure performed on clinical specimens is directly influenced by the quality of the specimens themselves. It is strongly recommended that laboratories employ methods discussed in the Manual of Clinical Microbiology¹⁷ for specimen collection, transport, and placement on primary isolation media.

Due to variations in inoculum concentrations prepared with McFarland standards, use of the BD BBL CrystalSpec Nephelometer, the BD PhoenixSpec Nephelometer, or the BD Phoenix AP instrument is required for adjusting the test inoculum prior to use in the BD Phoenix system.

Instructions for an optional purity check are provided at the end of this subsection.

Materials Required:

- BD Phoenix Panels
- BD Phoenix ID Broth or BD Phoenix Inoculum Broth
- BD Phoenix AST and/or BD Phoenix AST-S Broth
- BD Phoenix AST Indicator Solution and/or BD PhoenixAST-S Indicator Solution
- BD Phoenix Panel closures
- BD Phoenix Inoculation Station
- BD Phoenix Transport Caddy
- BD BBL CrystalSpec Nephelometer, the BD PhoenixSpec Nephelometer, or the BD Phoenix AP instrument
- 25 µL pipettor and tips

Materials Required but Not Provided:

- Gram Stain Reagents
- Sterile Cotton Swabs, Inoculation Loops or Needles
- Nonselective Culture Plated Media (see Section 4.3 Preparing Panels)
- Incubators
- Biohazard Disposable Container
- Markers, etc.
- Vortex mixer

NOTES

Exercise care in handling BD Phoenix panels. Handle panels by the sides only to avoid marking, smudging, or obscuring the bottom or top of the panel in any way.

Accession bar code labels affixed to a BD Phoenix panel:

Must not be of fluorescent material.

Should not cover any BD Phoenix panel reaction wells.

Should not cover the BD Phoenix sequence number (panel) barcode.

The procedure that follows describes all the steps in preparing a combination panel for both identification and susceptibility testing. If a combination panel is being used for only ID or only AST testing, note that certain steps are not applicable in the procedure.

General Panel Preparation

If the BD Phoenix AP instrument is being used, refer to the BD Phoenix AP Instrument User's Manual for panel preparation.

BD Phoenix Strep panels, BD Phoenix Yeast ID panels, BD Phoenix Emerge panels, BD Phoenix Inoculum Broth, and BD Phoenix MIC panels have separate instructions that appear after Step 19 of these instructions.

- 1 Confirm the Gram stain reaction of the isolate before proceeding with the inoculum preparation for use in the BD Phoenix M50 instrument. Once the Gram stain reaction is confirmed, select the appropriate BD Phoenix panel for inoculation.
- 2 Examine the pouch, and do not use the panel if the pouch is punctured or opened. Remove the panel from the pouch. Discard the desiccant. Do not use the panel if there is no desiccant or if the desiccant pouch is torn.

NOTE

Panels must be inoculated within 2 hours of being removed from the pouch.

- 3 Place the panel on the Inoculation Station; inoculation ports on top, and the pad on the bottom.

-
- 4 Label a BD Phoenix ID broth tube with the patient's specimen number. Using aseptic technique, pick colonies of the same morphology with the tip of a sterile cotton swab (do not use a polyester swab) or a wooden applicator stick from one of the recommended media.
 - 5 Suspend the colonies in the BD Phoenix ID broth (4.5 mL).
 - 6 Cap the tube and vortex for five seconds.
 - 7 Allow approximately ten seconds for air bubbles to surface. Tap the tube gently to aid in eliminating bubbles.
 - 8 Insert the tube into the BD BBL CrystalSpec or BD PhoenixSpec Nephelometer. Make sure the tube is inserted as far as it will go. (Refer to the BD BBL CrystalSpec or BD PhoenixSpec Nephelometer product insert for correct usage instructions.)
 - 9 If the inoculum density is set to 0.5 McFarland for the panel type being run, then a range of 0.50–0.60 is acceptable. If the inoculum density is set to 0.25 McFarland for the panel type being run, then a range of 0.20–0.30 is acceptable. If the density of organisms is low, colonies can be added from the isolate. Re-vortex the sample and reread to confirm that the correct McFarland has been achieved. If the density of organisms exceeds 0.6 McFarland, follow the steps below to dilute the broth. It is very important to accurately indicate the level of the liquid in the tube since this volume is needed to adequately fill the wells in the panel.
 - a Using a marker, mark the broth level in the over-inoculated BD Phoenix ID Broth tube.
 - b Using a sterile pipette, aseptically add fresh BD Phoenix ID Broth to the inoculum. Only BD Phoenix ID Broth may be used to dilute the inoculum.
 - c Vortex the tube and allow to sit for 10 seconds.
 - d Place the tube in the nephelometer and remeasure the turbidity of the suspension.
 - If the reading is greater than 0.6, repeat **Steps b-d**.
 - If the reading is 0.5 - 0.6, go to **Step e**.
 - e Using a sterile pipette, aseptically remove excess broth to the original level indicated by the mark on the tube created in Step a.

Remove excess broth to avoid overfilling the panel. Also, do not remove too much broth, as there may be insufficient broth to adequately fill the panel.
 - f Broth may now be used to inoculate the BD Phoenix AST Broth and/or the BD Phoenix Panel.

NOTES

- Yeast ID panels must be inoculated using a 2.00–2.40 McFarland inoculum density.
- Confirm current instrument settings for inoculum density before inoculating panels.
- See instructions below, ID Inoculum Density Flexibility, for information on using alternate densities.
- Only the BD PhoenixSpec Nephelometer and BD Phoenix AP instrument can be used to make inoculum densities of 0.25 McFarland
- Standardized bacterial suspension in ID Broth or Inoculum Broth must be used within 60 minutes of preparation.

- 10 If identification only is being performed, proceed to Step 15 and continue the procedure. If a BD Phoenix Emerge Panel is being inoculated, refer to the section below, BD Phoenix Emerge Panels.
- 11 Label a BD Phoenix AST broth tube (8.0 mL) with the patient's specimen number. Add one free falling drop of AST Indicator solution to the AST broth tube. Invert to mix. DO NOT VORTEX.

NOTES

- Allow AST Indicator Solution to warm to room temperature before dispensing into AST broth.
- The unused portion of the indicator should be returned to 2–8 °C as soon as possible. Do not store at room temperature for more than 2 hours. Opened bottles should be discarded after 14 days from initial opening.
- If volume other than one drop is added inadvertently, discard the tube and use a fresh tube of AST broth.
- After the addition of the Indicator to AST broth, the mixed solution can be stored in the dark, at room temperature, for as long as 8 hours.
- Tubes must be used within 2 hours after the addition of AST Indicator Solution if exposed to light.

- 12 If an inoculum density of 0.50–0.60 was used, transfer 25 µL of the bacterial suspension from the ID tube into the AST broth tube. If an inoculum density of 0.20–0.30 was used, transfer 50 µL (use 2 shots if utilizing a 25 µL pipettor) of the bacterial suspension from the ID tube into the AST broth tube.

NOTE

Panels must be inoculated within 30 minutes of the time that the AST broth inoculum is prepared.

- 13 Cap the AST tube and invert several times to mix.
- 14 Wait a few seconds for air bubbles to surface. Tap the tube gently to aid in eliminating bubbles.
- 15 Pour the ID tube inoculum into the fill port on ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel. If an AST (only) panel is being used, DO NOT inoculate the ID side of the panel. Retain the ID tube for an optional purity check (see below).
- 16 Pour the AST broth inoculum into the fill port on AST side of the panel (85-well side). Allow the fluid to traverse down the tracks before moving the panel.

-
- 17 Before placing panel closures check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present remove the droplet with absorbent material. The used absorbent material must be decontaminated before discarding.
 - 18 Snap on the panel closures. Make sure that the closures are fully seated. Use 2 closures regardless of panel type.
 - 19 Visually inspect panels to be sure each of the wells is full. Look at both sides of the panel. Make certain that the wells are not overfilled. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

See note above

- Panels must be loaded into the instrument within 30 minutes of inoculation.
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
- Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK

It is highly recommended that the purity of both ID and AST inocula be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24–48 h at 35 °C under appropriate conditions.

BD Phoenix Strep Panels

BD Phoenix Strep panels are for the identification and antimicrobial susceptibility testing of most *Streptococcus* species. Although *Streptococcus* species may be identified in the gram positive panels, antimicrobial susceptibility cannot be reported when using these panels. The BD Phoenix Strep panels, which must be used with BD Phoenix AST-S Broth and BD Phoenix AST-S Indicator Solution, provide the conditions required for rapid AST testing of most *Streptococcus* species.

- 1 Follow **Steps 1–9** (General Panel Preparation) to prepare the suspension of bacteria.
- 2 Add one drop of the BD Phoenix AST-S Indicator to each AST-S broth tube. Invert to mix. **DO NOT VORTEX.**

NOTES (See note above)

- **AST-S Broth and AST-S Indicator Solution are for use with the BD Phoenix Strep panels (SMIC/ID, SMIC) only. These reagents are not interchangeable with the AST Broth and AST Indicator Solution used with BD Phoenix Gram Positive and Gram Negative panels.**
- **Allow AST-S Indicator Solution to warm to room temperature before dispensing into AST-S broth.**
- **The unused portion of the indicator should be returned to 2–8 °C as soon as possible. Do not store at room temperature for more than two hours. Opened bottles should be discarded after 14 days from initial opening.**
- **If volume other than one drop is added inadvertently, discard the tube and use a fresh tube of AST-S broth.**
- **After the addition of the Indicator to AST-S broth, the mixed solution can be stored in the dark, at room temperature, for as long as 8 hours.**
- **Tubes must be used within two hours after the addition of AST-S Indicator Solution if exposed to light.**

- 3 Using a pipettor, transfer 25 µL of the standardized bacterial suspension from the ID tube into the AST-S broth tube.

NOTE

Panels must be inoculated within 30 minutes of the time that the AST-S broth inoculum is prepared.

- 4 Cap the AST-S tube and invert several times to mix.
- 5 Wait a few seconds for air bubbles to surface. Tap the tube gently to aid in eliminating bubbles.

-
- 6 Pour the ID tube inoculum into the fill port on ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel. If a BD Phoenix Strep MIC only panel is being used, DO NOT inoculate the ID side of the panel. Retain the ID tube for an optional purity check (see below).
 - 7 Pour the AST-S broth inoculum into the fill port on AST side of the panel (85-well side). Allow the fluid to traverse down the tracks before moving the panel.
 - 8 Before placing panel closures, check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present, remove the droplet with absorbant material. The used absorbant material must be discarded with biohazard waste.
 - 9 Snap on the panel closures. Make sure that the closures are fully seated.
 - 10 Visually inspect panels to be sure each of the wells is full. Look at both sides of the panel. Make certain that the wells are not overfilled. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

- Panels must be loaded into the instrument within 30 minutes of inoculation
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
- Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK


It is highly recommended that the purity of both ID and AST-S inocula be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24–48 h at 35 °C under appropriate conditions.

ID Inoculum Density Flexibility

An ID portion of a panel may be run in the opposite mode from what is configured by darkening well A-17 on the back of a panel before placing the panel in the instrument. This allows a panel to be run at an inoculum density of 0.20–0.30 even if a density of 0.5 is configured for that particular panel type. Likewise, a panel can be run at an inoculum density of 0.50–0.60 if a density of 0.25 is configured.

There is no way to alter the density setting during Panel Login. To use a panel in the opposite density mode, using a black permanent marker blacken the entire A-17 well as shown in the figure below.

Inoculum densities sent from the BD Phoenix AP instrument (via the BD EpiCenter system) cannot be changed by darkening well A-17. Inoculum densities sent from the BD Phoenix AP instrument are invalid if they are received after the completion of the first test cycle, or if the ID portion of the panel is disabled, or if the panel type does not support the inoculum density. If the instrument receives an invalid inoculum density for an ID panel that differs from the panel's determined inoculum density, the ID side of the panel fails, a Needs Attention code is triggered (Can Not Determine Organism ID), and a Special Message is triggered.

The following table outlines how to run the opposite mode from the density for which the instrument is configured.

Current Instrument Inoculum Density Configuration	Required Concentration (McFarland for Test Panel)	Amount of ID Inoculum to Add to AST Broth**	W-1A-17
0.50	0.25	50 µL	Blackened
0.25	0.50	25 µL	Blackened
** If also running AST			

BD Phoenix Yeast ID Panels

BD Phoenix Yeast ID panels are for the identification of most clinically relevant yeast and yeast-like species.

- 1 Follow **Steps 1–8** (General Panel Preparation) to prepare the suspension of yeast.
- 2 The inoculum density is set to 2.0 McFarland for Yeast ID panels, with a range of 2.00–2.40 as acceptable. If the density of organism is low or the density of organisms exceeds 2.40 McFarland, follow the steps discussed in Step 9 (General Panel Preparation) to obtain the correct McFarland density.
- 3 Pour the ID tube inoculum into the fill port on the ID side of the panel (51-well side). Allow the fluid to traverse down the tracks before moving the panel.
- 4 Before placing panel closures, check for residual droplets of inoculum on the edge of the fill ports. If a droplet is present remove the droplet with absorbant material. The used absorbant material must be discarded with biohazard waste.
- 5 Snap on the panel closures. Make sure that the closures are fully seated. Use 2 closures regardless of panel type.
- 6 Visually inspect panels to be sure each of the wells is full. If any of the wells are unfilled or overfilled, inoculate a new panel.

NOTES

See note above

- **Yeast ID panels must be inoculated using a 2.00–2.40 McFarland inoculum density.**
- **Standardized bacterial suspension in ID Broth or Inoculum Broth must be used within 60 minutes of preparation.**

NOTES

See note above

- Panels must be loaded into the instrument within 30 minutes of inoculation.
- Panels must be kept in the inoculation station after inoculation until the excess fluid has been completely absorbed by the pad.
- Panels should stay vertical in the caddy until loaded.
- Inoculated panels should be handled with care. Avoid knocking or jarring the panel.

NOTE

OPTIONAL PURITY CHECK

It is highly recommended that the purity of ID inoculum be checked by preparing a purity plate.

To perform a purity check, using a sterile loop, recover a small drop from the inoculum fluid tube either before or after inoculating the panel and inoculate an agar plate (any appropriate medium) for purity check. Discard inoculum fluid tube and cap in a biohazard disposal container. Incubate the plate for 24–48 h at 35 °C under appropriate conditions.

BD Phoenix Emerge Panels

BD Phoenix Emerge panels are designed to perform susceptibility testing on an expanded number of antimicrobial agents. To accomplish these susceptibilities, antimicrobial agents are present on both sides of the BD Phoenix panel. These panels do not have the ability to perform bacterial identification. Because of the design, the inoculation technique is unique and is outlined below. Two tubes of AST broth will be required.

- 1 Follow **Steps 1–9** (General Panel Preparation) to prepare the suspension of bacteria.
- 2 Add one drop of the BD Phoenix AST Indicator to each AST broth tube.
- 3 Transfer 25 μ L (50 μ L if low inoculum option is used) of the suspension to two BD Phoenix AST broth tubes. Cap and gently invert.
- 4 Using sterile technique, remove 3.5 mL of broth from one of the inoculated BD Phoenix AST broth tubes and discard in an appropriate container.
- 5 Pour the remaining 4.5 mL into the left side of the BD Phoenix Emerge panel. Pour the other BD Phoenix AST broth tube into the right side of the BD Phoenix Emerge panel.
- 6 Cap the panel and follow the normal panel login procedure.

BD Phoenix Inoculum Broth

BD Phoenix Inoculum Broth can be used to make the initial McFarland suspension of microorganisms when utilizing BD Phoenix MIC only panels (PMIC, NMIC, SMIC). The BD Phoenix Inoculum Broth is filled with 2.2 mL of BD Phoenix ID Broth and will reach the correct inoculum density by using fewer bacterial colonies.

- 1 Follow **Steps 1–9** (General Panel Preparation) to prepare the suspension of bacteria using the BD Phoenix Inoculum Broth rather than the BD Phoenix ID Broth.
- 2 For NMIC and PMIC panels, add one drop of the BD Phoenix AST Indicator to the BD Phoenix AST broth tube.
- 3 For SMIC panels, add one drop of the BD Phoenix AST-S Indicator to the BD Phoenix AST-S broth tube.
- 4 Transfer 25 μ L (50 μ L if low inoculum option is used) of the BD Phoenix Inoculum Broth suspension to the BD Phoenix AST or AST-S tube(s).
- 5 Follow **Steps 16–19** (General Panel Preparation). Follow the normal panel login procedure.

BD Phoenix MIC (NMIC, PMIC, SMIC) Panels

Inoculum can be prepared using either BD Phoenix ID Broth or BD Phoenix Inoculum (step 4, General Panel Preparation)

- 1 Follow **Steps 1–9** (General Panel Preparation) to prepare the suspension of bacteria.
- 2 For NMIC and PMIC panels, add one drop of the BD Phoenix AST Indicator to the BD Phoenix AST broth tube.

For SMIC panels, add one drop of the BD Phoenix AST-S Indicator to the BD Phoenix AST-S broth tube.
- 3 Transfer 25 μ L (50 μ L if low inoculum option is used with NMIC or PMIC panels) of the BD Phoenix ID Broth to the BD Phoenix AST or AST-S tube. Cap and gently invert. Retain the ID tube or Inoculum tube for an optional purity check (see above).
- 4 Follow **Steps 16–19** (General Panel Preparation). Follow the normal panel login procedure.

4.3.1 Quality Control

Refer to the panel package insert for information on ID/AST Quality Control. Note that QC panels cannot be marked as critical.

Quality Control testing is recommended for each lot of panels. The QC Lot Support feature can facilitate QC panel tracking and testing. If the QC Lot Support feature is enabled, then the panel lot number must be defined prior to logging in QC panels.

- 1 Inoculate a panel with one of the organisms listed in the package insert. All microbial cultures, including QC organisms, are potentially infectious and should be treated with universal precautions. QC organisms are prepared for panel inoculation as specified in Section 4.3 Preparing Panels, immediately preceding this section.
- 2 Log the panel in as a QC panel as follows:

NOTE

For most reliable results, it is recommended that the QC organisms be subcultured at least twice on two consecutive days onto TSA with 5% Sheep Blood agar before use in the BD Phoenix M50 instrument. Use only well-isolated colonies.

For Yeast ID QC, the use of Sabouraud Dextrose Agar as the subculture media is also acceptable.

- a Select **Panel Login**.
- b Select **QC**.
- c If desired, enter Tech ID.
- d If desired, in the Accession Number field, type in or scan an accession number.
- e In the Sequence Number field, type in or scan the panel's sequence number. If an invalid Sequence Number is entered, message W100 appears. Invalid Panel Sequence Number (see The field above does not display when BD EpiCenter communications is enabled.).
- f The Isolate Number field defaults to isolate number 1. Type in the isolate number, or select the + or the – to increase or decrease the number. Valid isolate numbers are 1 to 20.
- g For Yeast ID panels, the Media field appears to the right of the isolate number field. If a media type is not specified, a workflow error is generated when the panel is attempted to be saved. If a Yeast ID panel is not logged in before placing it in the instrument for testing, the panel aborts after the first reading because no media has been specified. A default media type can be configured (see 3.10.2 System Configuration Sub-Tab in Section 3.10 Configuration Tab), which appears when a Yeast ID panel sequence number is scanned during login. To select a different media, select the down arrow to display a drop down box listing all media types (abbreviations), sorted alphabetically. (When the media type is selected, the full name appears at the top right of the screen.) Highlight the desired media. Select **Save**.
- h In the Test Strain field, select the QC organism selection box from the drop down menu. Only the predefined, required test strains appear [if available], otherwise all of the test strains appear. If a Sequence number has not been entered, No Data Available is displayed. Highlight the desired QC organism. Select the organism.

- i In the Panel Lot Number field, type in or scan the panel's lot number. Lot numbers must be seven numeric digits. Press tab to advance to the next field. (If the QC Lot Support feature is enabled, this field is completed automatically when the Sequence Number barcode is scanned.)
- j In the Expiration Date field, select the calendar icon to enter the expiration month, day, or year. (If the QC Lot Support feature is enabled, this field is completed automatically when the Sequence Number barcode is scanned.)
- k If desired, in the ID Broth Lot Number field, type in or scan the broth lot number. Lot numbers must be seven numeric digits.
- l If desired, in the Expiration Date field, press the calendar icon to select the expiration month, day, or year
- m If desired, in the AST Broth Lot Number field, type in or scan the broth lot number. Lot numbers must be seven numeric digits.
- n If desired, in the Expiration Date field, press the calendar to select the expiration month, month, or year.
- o If desired, in the Indicator Lot Number field, type in or scan the broth lot number. Lot numbers must be seven numeric digits.
- p If desired, in the first Expiration Date field, press the calendar to select the expiration month, day, or year.
- q Select **Save** to save the information.

Place the panel in the instrument (refer to Section 3.3.2 Inserting Panels in the Instrument).

Figure 4-1 – QC Panel Login Screen

The field above does not display when BD EpiCenter communications is enabled.

4.4 Automatic Association of Panels

4.4.1 Overview

The purpose of Automatic Association is to associate ID results from one panel (an ID or Combination panel) to the AST panels that are related to it and that lack an ID. Panels are related by virtue of having the same accession and isolate number. (Thus panels with the same accession number and different isolate numbers are NOT related, and Auto Association does not occur between such panels.)

Auto Association often involves the initial stage of Organism ID Conflict Checking. The circumstances where Conflict Checking occurs are described below. Conflict Checking is the process where the system verifies that the panel which has triggered the process does not have a related panel with a Final ID that is different from its own. If there is a conflict, Auto Association does not occur between related panels.

Auto Association can help eliminate unnecessary work such as the need to manually enter an organism ID for related panels. In addition, the Conflict Checking function helps ensure that the same organism ID is made for all the panels related to a patient specimen.

Note that Auto Association does not set the critical panel attribute for related panels.

4.4.2 Panel Types

The panel types whose ID information is associated to related panels are: ID panels, Combination panels (with both ID and AST sides enabled), and Combination panels with only the ID side enabled. Auto Association uses the Final ID.

The related panels whose IDs are set by Auto Association are: AST panels, and Combination panels with only the AST side enabled.

4.4.3 When Auto Association Is Not Performed

Auto Association is NOT performed in the following circumstances:

- When an organism ID conflict exists within an accession/isolate
- To a related AST panel that has been finalized or deleted, or that already has a Final ID
- Between related ID-type panels
- Between related AST-type panels
- Between QC and regular panels
- Between related QC panels
- If the Instrument ID field contains a tie/triplet, this information is not associated to any related AST panels, however the Final ID information IS auto associated

4.4.4 Organism ID Conflict Checking

Organism ID Conflict checking is performed using the Final ID for panels within the same accession/isolate (i.e., related panels) that are not finalized (and are not QC panels). If both panels have a Final ID (both are not blank) and they do not match exactly, an Organism ID Conflict is considered to exist. When an Organism ID Conflict exists, no Auto Association is performed. All the panels involved in the conflict have a Needs Attention flag set indicating that there is an organism conflict. Once the conflict is resolved (and provided that all other conditions for Auto Association are met), then Auto Association does occur.

Note that Organism ID Conflict checking is not performed on QC panels.

4.4.5 Auto Association and Related Actions

The trigger events that invoke Auto Association include: an ID or Combination panel status transitioning to Removable or accession/isolate information being modified. If there are no unresolved tie/ triplet conditions and no Organism ID Conflicts exist:

- 1 The instrument searches for the first unfinalized related ID or Combination panel with a Final ID. If a panel with these criteria is found, the instrument performs Auto Association.
- 2 If a panel with the above criteria is not found, the instrument then searches for the first finalized related ID or Combination panel. If a finalized panel is found, the instrument performs Auto Association.
- 3 If there are related ID panels with no Final ID when an AST panel invokes Auto Association (and no Organism ID Conflict condition exists), the instrument sets a Missing Organism ID Needs Attention reason for the AST panel if its status is Removable.
- 4 Auto Association is performed for any related AST panels that do not have information in their Final ID fields if the AST panel's status is either ongoing or Removable (Auto Association is not performed on AST panels with a status of Pending).
- 5 Auto Association is the action of modifying an AST panel's Final ID field (when this field is blank) to contain the same information as a related ID panel's information.
- 6 The instrument ID information is only associated provided it is a single Organism ID and not a Tie/Triplet condition.
- 7 A Combination panel with both sides enabled uses its ID result to perform interpretations on the AST side of the panel.
- 8 When Auto Association is performed successfully, any existing Missing Organism ID Needs Attention reasons are automatically cleared from the AST panels.

4.4.6 Typical Auto Association Examples

Instrument ID Is Associated

ID and AST panels are logged with the same Accession/Isolate number. (The AST panel is not logged in with an Organism ID.) They are both placed into the instrument at the same time. The ID panel finishes testing. The panels look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A
AST		

When the ID panel has completed testing, it looks for related panels. It sees the AST panel with no Final ID. The ID panel checks for Organism ID Conflicts. There is no Organism ID Conflict condition, so the ID panel associates its Instrument ID to the AST panel record's Instrument ID field (regardless of the AST panel's status).

The panels look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A
AST	Org A	Org A

Instrument ID and Final ID Are Associated

The ID panel is logged in with an Accession/Isolate number. It is placed into the instrument to test. It completes testing and looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org A

The user changes the Final ID for this panel. It now looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org B

The user then logs in an AST panel with the same Accession/Isolate number as the ID panel. (The AST panel is not logged in with an Organism ID.) The AST panel completes testing and looks for related panels. It sees the ID panel and checks to make sure that there is no Organism ID Conflict within the Isolate. There is no Organism ID Conflict, so the AST panel associates the ID panel's ID to itself. The panels now look like the following:

Panel	Instrument ID	Final ID
ID	Org A	Org B
AST	Org A	Org B

Final ID Is Associated

An ID panel is logged in with an Accession/Isolate number. The panel is placed in the instrument to test. The ID panel finishes testing first with a tie /triplet condition. When this situation occurs, the panel looks like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	

When the tie/triplet is resolved (by the user selecting one of the organisms listed or another organism), the ID panel will look like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	Org A

The user then logs in an AST panel with the same Accession/Isolate number as the ID panel. (The AST panel is not logged in with an Organism ID.) The AST panel completes testing and looks for related panels. It sees the ID panel and checks to make sure that there is no Organism ID Conflict within the Isolate. There is no Organism ID Conflict, so the AST panel associates the ID panel's Final ID to itself. The panels now look like the following:

Panel	Instrument ID	Final ID
ID	Org A	
	Org B	
	Org C	Org A
AST		Org A

4.4.7 BDxpert Triggered Rules Screen

By default, the BDxpert Rules are enabled and set to trigger automatically when test results are processed. However, there is an option in the BDxpert Rule Configuration to disable some rules (which means they would not be applied to results) or to set BDxpert rules to trigger manually. Any rules set to trigger manually must be accepted or rejected by the user before any additional rule processing or Final SIR determination occurs.

To access BDxpert Triggered Rules:

- 1 Select **BDxpert Rules** from the Results screen.
- 2 The first pending rule is highlighted.
- 3 The effect of the rule can be viewed on Final SIR values in the AST Results screen.
- 4 Select **Accept** to accept the pending rule, **Reject** to reject the rule, and **Re-Run** to cause all interpretation rules to be reapplied to the raw data.

4.4.8 Resistance Markers

If Resistance Markers have triggered on a particular panel, they are listed with the BDXpert rules.


The following information is shown:

- The BDXpert rule number that triggered the Resistance Marker
- The Resistance Marker code (abbreviation)
- The Resistance Marker Name
- The BDXpert rule Description: Accept, Reject, Re-execute or Exit

5 – Maintenance

5.1 General

The BD Phoenix M50 instrument requires little maintenance from the user to provide reliable performance. Daily activities include checking the instrument temperature and printer paper supply. On a weekly basis check the operation of the station status indicators (LEDs), the audible alarm, and the alert indicator. Routine preventive maintenance consists of a monthly check of the ambient air filter. All other procedures are on an as needed basis. Any maintenance or repair not described in this section should be performed by BD personnel only.

No yearly preventative maintenance is required to be performed by BD authorized service personnel.

5.2 Routine Maintenance

Time Frame	Procedure
Daily Record results on Daily Instrument Report	• Check paper supply. • Record temperature on the Status screen. The temperature should be $35\text{ }^{\circ}\text{C} \pm 1.5\text{ }^{\circ}\text{C}$. • Record temperature standard panel. The temperature standard panel can be brought into view by selecting one of the instrument LED check functions on the Maintenance tab. • The temperature should be $35\text{ }^{\circ}\text{C} \pm 1.5\text{ }^{\circ}\text{C}$. • Perform Daily Verification on BD PhoenixSpec. • See Nephelometer Calibration in the BD PhoenixSpec package Insert. The results are recorded on the Daily Instrument Report.

Time Frame	Procedure
Weekly	• Test Internal Green LEDs • Test Internal Red LEDs • Test Internal Amber LEDs • Extinguish all LEDs • Test External System LED • Test Alarm Perform the tasks listed under the Maintenance Hardware category. See below for description of task.
Every 3 Months	Perform calibration of BD PhoenixSpec. • See Nephelometer Calibration in BD PhoenixSpec package Insert. Clean panel stations if they become contaminated by a leaking panel. The priority in this situation is to first limit the extent of the contamination and then to decontaminate the panel location(s) and other accessible instrument areas receiving the spill. • The contamination procedure applies only to accessible areas. If the spill extends into regions of the carousel that are not accessible, contact your local BD Service representative.
As Needed	• The solution recommended to clean the affected surfaces should be at least a 10 percent household bleach solution. All surfaces must be thoroughly washed with the freshly prepared bleach solution, so that the surfaces are glistening wet. If the extent of the contamination is uncertain, thoroughly wash the exposed portions of the carousel and cabinet with the freshly prepared bleach solution. See Section 5.2.5 Cleaning / Decontamination for a detailed cleaning procedure.

Time Frame	Procedure
	Replace the Barcode Scanner 1 Locate the cable connecting the barcode scanner to the instrument. 2 Unplug the cable. 3 Plug in the new scanner cable. 4 Verify the proper operation of the new scanner by scanning a panel sequence number at the Login screen. 5 View that the sequence number is correct. • Select Cancel to exit the screen without saving the panel.

5.2.1 Cleaning and Checking the Air Filter

The filter should be checked monthly and cleaned/replaced if needed (see below, Cleaning the Air Filter).

If the instrument's environment is especially dusty, the air intake filter should be checked more frequently and cleaned or replaced if needed. The filter must remain clean and unobstructed; restricted airflow from a dirty filter may cause the instrument interior to reach excessive temperatures, which can affect results and possibly cause hardware malfunctions or failures. The filter can be cleaned and reused.

The instrument's filter is located inside the panel accesses area on the lower right. The filter can be removed without tools.

5.2.2 Removing the Air Filter


Figure 5-1 – Air Filter Location

- 1 Press Panels in or Panels out button to access panel entry area.
- 2 The filter is located on the lower right side of the opening.
- 3 To remove the filter, grasp black tab and slide towards the left side of the panel access area.
- 4 Remove the old filter and clean and dry it before replacing in the instrument, or place a new filter in the housing while the old one is cleaned and dried.
- 5 To insert a clean filter, grasp the black tab on the side of the filter and slide to the right until filter flange is flush with right side of instrument access area.
- 6 Close the instrument door.

5.2.3 Cleaning the Air Filter

- 1 Wash the dirty filter in a bactericidal disinfectant.
- 2 Place the filter on a paper towel and dry it thoroughly (if it is going to be reused immediately).
- 3 To save time, replace the dirty filter with a spare clean filter. Wash, dry, and set aside the removed dirty filter for the next filter replacement.


Figure 5-2 – Air Filter

5.2.4 Daily Instrument Report

The Daily Instrument Report can automatically print at a user defined time or it can be printed upon request. (see Figure 5-3 – Daily Instrument Report).

For instructions on configuring the report to automatically print at a certain time, see the Instrument Configuration section of this manual.

The Daily Instrument Report has three sections:

Section	Description
Top	Lists information that was configured when the system was set-up: • Facility name and address • Date and Time the report was printed • Instrument Number • Serial Number
Middle	Lists information about the internal operation of the Phoenix as well as the Normalizer Panels. • This section should be reviewed on a daily basis to make sure that these operations have passed testing.
Bottom	Lists the maintenance items that need to be checked and recorded on a Daily and Weekly basis.

Daily Instrument Report

The Diagnostic Switch Is ON

7/8/2018 3:20 PM
1.0.65.0 / V5.91A (x-US)

System Number56

AB Serial NumberPFVV0015

Instrument TemperaturePass

Carousel Rotational TestPass

Power Supply Check?

Normalizer Panels

Tier	Sequence Number	Status	Expiration Status
A	429932154527	Pass	7/9/2018
B	429932154528	Pass	7/9/2018

Daily	Pass	Fail	Tech ID
Record Instrument Temperature From Main Screen Range: (33.5 - 38.5)			
Record Standard Panel Temperature Range: (33.5 - 38.6)			
Check Printer Paper Supply			
Weekly			
Instrument Internal Green LEDs			
Instrument Internal Red LEDs			
Instrument Internal Amber LEDs			
Instrument System Alert Indicator			
PC Audible Alarm			
Comments			

This report may contain PHI and/or PII, handle appropriately.

Page 1 of 1

Figure 5-3 – Daily Instrument Report

5.2.5 Cleaning / Decontamination

A situation requiring biological decontamination of one or more panel locations can occur if a panel should leak while in the instrument. The priority in this situation is to first limit the extent of the contamination and then to decontaminate the panel location(s) and other accessible instrument areas receiving the spill. If the spill extends into regions of the carousel not accessible for topical decontamination, contact your local BD representative for further instructions.

To Decontaminate Carousel Panel Locations


The solution recommended to clean the affected surfaces should be at least a 10 percent household bleach solution. All surfaces must be thoroughly washed with the freshly prepared bleach solution, so that the surfaces are glistening wet. If the extent of the contamination is uncertain, thoroughly wash the exposed portions of the carousel and cabinet with the freshly prepared bleach solution.

Required Materials:

- 10% bleach solution
- Personal protection equipment, including gloves, gown, eye protection (e.g., face shield, goggles, etc.)
- Gauze pads or paper towels
- Tap water

Cleaning Procedures:

- 1 Wear gloves and a gown, completely covering any body surfaces that could possibly come into contact with the affected instrument surfaces.
- 2 Turn off power to the instrument. Unplug the instrument power cord before proceeding.
- 3 Completely absorb the contaminated spill (gauze pads are most effective).
- 4 Apply the bleach solution to the affected surfaces, so that the surfaces are glistening wet. Let stand for approximately 15 minutes.
- 5 Absorb the applied solution with gauze pads or paper towels.
- 6 Dampen a clean cloth with water. Wipe down the decontaminated surfaces.
- 7 Thoroughly dry all wet surfaces.
- 8 Discard all cleanup materials as biohazardous waste.

5.3 Module Replacement

5.3.1 General

The BD Phoenix M50 instrument has been designed and tested for trouble-free performance. In the event of malfunction, contact BD for service under existing contract terms or warranty. Only the external barcode scanner is user-replaceable.

Replacement modules may be exchanged. Only replacement parts supplied by BD should be used in the procedures described in this section.

5.3.2 Thermometer Removal

If the fluid in the thermometer of the temperature standard panel has separated, follow the procedure below to remove the thermometer and replace it.

To Remove/Replace the Thermometer

- 1 Remove the temperature standard panel from the BD Phoenix M50 instrument.
- 2 Using the small access hole at the bottom of the panel, gently push the thermometer upward through the large slotted opening at the top of the temperature panel.
- 3 Manually pull the thermometer completely out of the temperature standard panel through the slotted opening at the top of the temperature panel.
- 4 Reunite the separated fluid column per the instructions below.
- 5 Install the thermometer in the reverse order described above.

5.3.3 Reuniting Separated Liquid in the Thermometer

If the fluid in the thermometer of the temperature standard panel has separated, follow the procedure below to reunite the liquid.

To Reunite Separated Liquid

- 1 Should there be a separation in the capillary or in the expansion chamber at the top of the thermometer, heat the bulb of the thermometer in a hot liquid which exceeds the range of the thermometer until the separation and main liquid column enter the expansion chamber and unite with each other.
- 2 Quickly remove the thermometer from the liquid, so that the liquid does not completely fill the expansion chamber which could possibly harm the thermometer.
- 3 Check the thermometer against a certified, traceable thermometer or an ice bath to assure the thermometer is reading correctly.

6 – System Alert Messages

6.1 Error / Alert Messages

CAUTION

When the system displays alerts and errors, immediately respond to the condition.

When the system encounters an alert or error condition, the error code (EXX or WXXX, where XX or XXX is a number) is either displayed on the screen or written into the system alert list. The error code is an abbreviation for the conditions described in the listing below.

Different types of alerts and errors behave in different ways. There are three basic types of alert conditions:

- Self-Clearing – these alerts are removed from the System Alert Screen after they have been displayed and the screen has been exited. No other intervention is required.
- Persistent – these alerts remain in the System Alert Screen (even after displaying the alert code) until the instrument determines that the error causing the alert has been corrected. Correction can be accomplished with or without user intervention, depending on the alert.
- Auto Clear – these alerts are removed from the System Alert Screen as soon as the condition causing the error has been cleared. There is no requirement to display the alert code on the System Alerts List Screen. If the alert is still present when reviewing the Alert list, it will be displayed.

Each error code called out in this chapter will contain one of these three types of alerts and errors.

W error codes are displayed on the screen when they occur. (They also cause the Workflow Error tone to sound [sequence of short high beep and short low beep repeated four times].) These are activity (or workflow) types of errors. In most cases, this means that some action that was performed was not what the system expected, but the correct action can usually be performed, as recommended below, without exiting the current operation. These activity errors are flagged by a Workflow Error (see Figure 6-1 – Workflow Error).


Figure 6-1 – Workflow Error

W error codes are grouped as follows and full descriptions can be found in Section 6.3 Workflow Alerts (W error codes):

- W1XX Problem with the Sequence Number
- W2XX Problem with the Accession, Isolate Number, or Media type
- W3XX Problem with a QC panel only
- W4XX Action to get into the instrument is not allowed
- W5XX Configuration/Maintenance screen activity not allowed (as long as it is not an instrument action)
- W6XX Screen activity is not allowed
- W7XX LIS errors
- W8XX QC Lot Support errors

6.2 System Alerts (E error codes)

System alerts (each shown in separate tables below), comprise of all E type error codes and are reported in the system alert list. These errors cause the Alert tone (medium beep on for one second, off for 3 seconds, repeating) to sound (if it is enabled). Also the System Alert icon appears on the Status screen and the Alert indicator is illuminated. The errors must be reviewed to clear the system alert condition. The system alert list can be viewed from the Status screen by touching the alert icon on the upper left corner of the display.

The E error codes are listed in numerical order. Error sub-codes are 8-digit numeric codes that appear below the EXX readout in the system alert list. The sub-codes indicate specific conditions detected and many are listed in the associated alert tables below

CAUTION

If any error sub-codes other than those listed here appear, note the sub-code and contact BD for assistance. If the recommended corrective actions do not solve the problem, contact BD.

E01 Incubator Temperature			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	Average temperature too high (> 36.5 °C)	For all sub-codes, causes include: room temperature is not within recommended range, or other environmental specification is not being met (such as instrument sitting in direct sunlight or too close to HVAC air register). Make sure environmental specifications are met (see Environmental Requirements in Section 2.3 Instrument Specifications) Clean and/or repair air filters to permit fresh air intake.
00000002		Average temperature too low (< 33.5 °C)	
00000004		Communication with the incubator temperature sensor is lost	
00000010		Average temperature too high for more than one hour (> 36.5 °C)	
00000020		Average temperature too low for more than one hour (< 33.5 °C)	
00000040		Absolute high temperature (> 38.5 °C)	
00000080		Control and QC temperatures disagree	
00000100		QC temperature is out of range	

Temperature readings are taken every 5 seconds for 10 minutes and these readings are averaged. Note that if the instrument temperature reaches 38.5 °C, the instrument disables the heater.

E05 Carousel Alert			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	Carousel RPM is below specification.	The instrument reports these alerts if something is impeding the motion of the carousel. This can be attributed to one of the following: the reported RPM is below spec, certain flag readings are incorrect, or the carousel is jammed or stalled. It is detected during any carousel rotation. Any inventory scan or panel test in progress is aborted and the instrument ignores any data received from the test or scan. Open the door, look for and remove any obstructions such as a panel that is ajar or a panel closure that is not seated. Do not manually rotate the carousel. Close the door. If message reoccurs, contact your local BD representative.
00000080		Carousel is jammed (no panel flags detected, drum not moving)	

Tier Alert			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	This is the general alert condition for problems related to specific tiers including: incorrect flag readings, optical errors, timeout conditions, and normalizer problems.	Tier alerts are detected during a test cycle or inventory scan. The instrument ignores any data received from the test or scan for the tier that has sustained the error. Data received from good tiers is retained. All locations in the bad tier are automatically blocked. Some conditions reported for the E05 alert are also detected by the tier micro and reported with the E06 alert. They include missing/extra panel flags, home flag, etc. Follow screen instructions to correct the alert.
00000002			
00000004			
00000008			
00000010			
00000040			
00000080			
00000100			
00000200			
00000400			
00004000			
00020000			
00040000			
00080000			
00100000			
00200000			
00400000			

* The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E09 Test Aborted Sub-code			
Sub-code	Alert Type	Possible Causes	Corrective Actions
20000000	Self-clearing	Panel testing has not occurred for more than one hour because the instrument was off, the door was open for more than an hour, the system clock was set ahead one hour or more or the test cycle had not occurred in more than an hour.	All ongoing panels are set to Needs Attention, and their status is set to Removable. Some results may be incomplete, and all affected panel results should be reviewed.

E10 System Database Corruption			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000002	Self-clearing	System Parameters Database Corrupted	Check your settings in the Configuration screens (see Section 3.10.2 System Configuration Sub-Tab) and reset them to your preferences. Save data to a USB key and contact your local BD representative.
00000008		BDXpert Rules Database Corrupted	Check your settings for BDXpert Rule Configuration (see Section 3.4.5 BDXpert Rules Sub-Tab) and reset them to your preferences. Save data to a USB key and contact your local BD representative.
00000020		Custom Breakpoint Database Corrupted	Check your breakpoint settings for Custom Interpretation Rule Set Configuration (see Section 3.10.6 Custom Interpretation Rule Set Sub-Tab) and reset them to your preferences. Save data to a USB key and contact your local BD representative.
00000040		User Codes Database Corrupted	Check your settings in the Organism and Antimicrobial Configuration screens see Sections 3.10.4 Organism Configuration Sub-Tab and 3.10.5 Antimicrobial Configuration Sub-Tab) and reset them to your preferences. Save data to a USB key and contact your local BD representative.
00000001 00000010 00000080 00001000 00002000 00004000 00008000 00010000 00020000 00040000 00080000 00100000 00200000 00400000		Database corruption: Panel, Panel Lot, Alert / Eventlog. Corruption of instrument configuration record, instrument history, system parameters, light source, etc.	Save data to a USB key and contact your local BD representative.

E11 Printer Error

Sub-code	Alert Type	Possible Causes	Corrective Actions
20000000	Self-clearing	Paper jam or power condition	Check printer paper (jammed our out), cable connection, power on, and/or online indicator.

E13 Power Failure

Sub-code	Alert Type	Possible Causes	Corrective Actions
0x20000000	Self-clearing	Power removed from instrument	Message is informational. If multiple power failures have occurred, only the latest one is reported in the alert list. Note the power failure and restore times in your instrument log. Note that power fail events are not recognized until the instrument user interface has successfully loaded.

E14 CCD Underrun*

Sub-code	Alert Type	Possible Causes	Corrective Actions
00000000	Persistent	Scanning of a panel stopped prematurely during a test cycle	The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked. To clear, re-boot the instrument.

*The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier

E18 Normalizer Row Averages*			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001 00000002 00000004 00000008 00000010 00000020 00000040 00000080	Self-clearing	The tier had a problem with normalizer panel data averages	A source adjustment can only be performed when all panels in an instrument are no longer ongoing. Normalizer errors are detected during a test cycle. Test data received may be discarded or retained depending on the error sub-code. All available stations in the bad tier are automatically blocked; ongoing stations become blocked as testing completes or panels are removed. Do not enter any new panels or move an ongoing panels into any tier reporting this error. The instrument may be able to correct this error via the light source adjustment process.

* The first line of the error represents Tier A, the second line represents Tier B, the third line represents Tier C, and the fourth line represents Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E20 Barcode Scanner Not Communicating			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001 00000002 00000004 00000008 00000010	Self-clearing	Tier barcode scanner (or handheld scanner) is not communicating and all the tier stations are automatically blocked	Reboot the instrument.

E21 Level 2 Rotor Step Warning*

Sub-code	Alert Type	Possible Causes	Corrective Actions
00000000	Persistent	The instrument has sensed that rotor rotation was not ideal but still acceptable	Check for a panel protruding from its carrier, an improperly seated closure, or user applied label peeling off. If no obvious visible cause for the error exists, contact your local BD representative.

* The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E22 Level 3 Rotor Step Error*

Sub-code	Alert Type	Possible Causes	Corrective Actions
00000000	Persistent	The instrument has sensed that rotor rotation was out of specification	The instrument ignores any data received from the test for the panel/station that has sustained the error. Data received from good panels is retained, unless the station sustaining the error was the normalizer station. All error stations are automatically blocked. If the normalizer has this error, the whole tier is blocked. Check for a panel protruding from its carrier, an improperly seated closure, or user applied label peeling off. If no obvious visible cause for the error exists, contact your local BD representative.

* The first line of the error represents stations with errors in Tier A, the second line represents stations with errors in Tier B, the third line represents stations with errors in Tier C, and the fourth line represents stations with errors in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E30 Normalizer Expiration Alert*			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Self-clearing	The tier's normalizer panel expiration date is between 60 and 30 days away	The instrument issues a weekly alert beginning when Normalizer panel expiration is 60 days away, which progresses to a daily alert when expiration is 30 days away. Schedule Normalizer panel replacement for the affected tiers before they expire (expiration date is shown on Daily Instrument report).
00000002		The tier's normalizer panel expiration date is less than 30 days away but has not expired	Contact your local BD representative.

*The first line of the error represents Normalizers in Tier A, the second line represents Normalizers in Tier B, the third line represents Normalizers in Tier C, and the fourth line represents Normalizers in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E31 Normalizer Expired*			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	Normalizer panel expiration date has passed. Stations in all affected tiers are blocked	Schedule normalizer replacement immediately to resume testing.

* The first line of the error represents Normalizers in Tier A, the second line represents Normalizers in Tier B, the third line represents Normalizers in Tier C, and the fourth line represents Normalizers in Tier D. A code of 00000000 indicates that there is no error condition for that tier.

E44 LIS Fatal Operating System Error

Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	The instrument is not able to send data to the LIS due to a fatal software exception that occurs within the LIS library at the instrument	Start LIS system. Look for obvious source of problem such as disconnected cable. If no obvious source exists, contact your local BD representative.
00000008		A fatal operating error was detected by the LIS_IM	Contact your local BD representative.

E50 Internal Software Error

Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Self-clearing	System encountered a software general protection error.	Save data to USB key (see Section 3.9.1 Maintenance Hardware Functions) and contact your local BD representative.
00000002		Internal software error	

E51 Duplicate Panel Barcode Detected			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000000	Self-clearing	Duplicate barcodes have been detected.	Carefully examine the reported stations for duplicate barcodes. All panels with duplicate barcodes have been aborted. If duplicates are found, then remove the panels from your system and contact your local BD representative immediately. If no duplicates are found, then the alert may be a result of moving an existing panel between instruments. Avoid moving panels across instruments as they will be aborted. If the alert reoccurs, contact your local BD representative.

E52 Instrument Communications Application Layer (ICAL) Alert			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Self-clearing	Instrument is not registered	Configure the instrument in the system.
00000002		Instrument history purged because instrument was in isolation mode too long.	Reconnect the instrument and Tablet PC or contact your local BD representative.
00000004		Extra Device Alert: three or more instruments are plugged into the application.	Remove the instrument(s) that are not configured in the system.

E98 PC Alert			
Sub-code	Alert Type	Possible Causes	Corrective Actions
00000001	Persistent	The IDS layer timed out	If the condition reoccurs, contact your local BD representative.
00000002		The IDS system received an invalid parameter	Contact your local BD representative.
00000004		The IDS system rejected a message	
00000008		An invalid alert sub-code was detected	

6.3 Workflow Alerts (W error codes)

Field level alerts are displayed when data, entered into a specific field, is invalid. These errors (such as attempting to enter an invalid Sequence Number in the Panel Login screen) cause workflow error messages to appear on the displayed screen. They do not put the system into an alert condition. These errors can frequently be cleared by simply performing the activity correctly (such as entering a valid Sequence Number). When the error is corrected, the field level alert is removed from the display. All other workflow alerts are displayed in a dialog with two or more buttons.

W100 Invalid Panel Sequence Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	The panel sequence number typed or scanned does not meet the required number format.	Panel sequence number barcode is located at the top of the reaction side of the panel. Scan or type in the correct panel sequence number. Check panel carton for panel update barcodes if they are present.

W101 Missing Panel Sequence Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	Panel sequence number field empty when save is pressed.	Type in or scan the correct panel sequence number before attempting the operation again.

W102 Unknown Panel Sequence Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save or select a panel sequence number that is not in the BD Phoenix database.	Verify that the correct panel sequence number was entered. If error reoccurs, the panel may need to be logged in if it is a new panel. Older, completed, finalized panels eventually age out of the database.

W103 Invalid Panel Type for Region		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to log in a panel barcode that is not valid for the region.	Verify that the correct panel barcode was entered.

W200 Invalid Accession Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made save, find, or print information for an invalid accession number.	Enter a valid accession number, up to 20 characters excluding : * ? [] ! #

W201 Missing Accession Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An operation (find, save, print) was attempted and the value in the Accession Number field was invalid or blank. This can include: a record with a valid panel Sequence Number and an Isolate Number greater than 1 with no Accession Number; an orphan panel with just an Isolate Number, or trying to change a saved record to a blank Accession Number.	Type in or scan the correct accession number before attempting the operation again.

W202 Unknown Accession Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	At attempt was made to find data or print a report for an accession number that is not located in the BD Phoenix database.	Verify that the correct accession number was entered. If error recurs, the record may need to be logged in if it is new. Older, completed, finalized records eventually age out of the database.

W203 Missing Isolate Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a record when an accession number is present without an isolate number.	To save a record, if you enter an accession number, you must also enter an isolate number

W204 Missing Media Type		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a Yeast ID panel during Panel Login and no Media is specified.	To save a record, if you enter a Yeast ID panel sequence number, you must select a media type.

W205 Panel Not Found		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to find a Sequence and Accession combination that is not in the BD Phoenix database.	The Sequence and Accession combination does not exist in the database.

W300 Missing Test Strain for a QC Panel		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to save a QC panel without an organism ID (Test Strain).	You must select an organism (Test Strain) to save a QC panel. Advance to the Test Strain field to drop down a box listing the available test strains. Highlight the desired organism. Select the highlighted organism.

W301 Missing Tech ID for a QC Panel		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a QC panel without a Tech ID.	You must enter a Tech ID to save a QC panel. Advance to the Tech ID field and enter a Tech ID, up to 3 alphanumeric characters excluding: * ? [] ! # .

W302 Invalid Tech ID for a QC Panel		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a QC panel with an invalid Tech ID.	Enter up to 3 alphanumeric characters excluding: * ? [] ! # for the Tech ID.

W303 Invalid Lot Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a record or print a report for an invalid panel lot number, ID broth lot number, AST broth lot number, or indicator lot number.	Enter the correct lot number. Lot numbers can be up to 7 characters. The lot number is shown on the item carton.

W304 Missing Lot Number		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save or print a report for a QC panel without a panel lot number, or an expiration date has been entered for one of the optional lot number fields (ID Broth, AST Broth, Indicator) with no corresponding lot number.	You must enter a panel lot number to save a QC panel. Advance to the Panel Lot # field and enter a lot number, up to 7 digits. The lot number is shown on the panel carton.

W305 ID or AST Must be Enabled		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to save a panel or QC panel, and both the ID and AST fields are unabled.	The ID/AST fields are enabled or disabled based on the type of panel, according to the panel sequence number. At least one field must be enabled for the panel record to be saved and for testing to occur.

W400 Door Already Open		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	A request was made to perform a task that requires the door to be opened when the door is already open.	Select OK.

W401 Can Not Perform Panel Locate

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	The panel you are trying to locate is not in the instrument, or you cannot access the instrument using the Normalizer Panel Replacement function (e.g., because there are ongoing panels).	Verify that the correct panel sequence number is entered or that the instrument can be accessed (for Normalizer Panel Replacement activity, ongoing panels render the instrument inaccessible for this activity).

W500 Upgrade Error

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	The USB key does not contain a readable, same, or newer version of instrument software. Causes include: older version of instrument software; USB key was removed before the update completed, or the USB key contains a corrupted or missing file.	Attempt the update operation again. If error reoccurs, contact your local BD representative.

W502 Removable Media Error

Message Type	Alert Cause	Alert Corrective Action
Retry Cancel Print Screen	The USB key is not in the drive, full, write protected or it was removed before completion of the task.	Attempt the update operation again. If error reoccurs, contact your local BD representative.

W503 Duplicate Rule		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	The data contained in the Antimicrobial, Org(anism) Group, and Organism fields is the same for the rule being saved as a different rule already defined in the Rule Set.	One of the parameters (Antimicrobial, Organism Group, or Organism) must be unique for the record to be saved.

W506 Instrument Not Idle		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt to execute an upgrade failed because the instrument was not idle.	Wait until instrument is idle (or correct carousel jam) before performing this activity.

W507 Invalid Password		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An invalid (or no) password was entered to access a Configuration screen.	Enter the correct password and perform the action again.

W511 Battery Mode

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt to install or restore data while the tablet is unplugged and working in battery mode.	Retry the install or restore data when battery mode is not in progress.

W600 Panel Can Not Be Deleted

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to delete a panel that is still physically in the instrument or a non-pending QC panel that is not physically in the instrument (if QC Lot Support is enabled).	Remove the panel from the instrument. You can use the find panel tab to make the carousel present the panel when the door is opened. If you know you have already removed the panel, then close the instrument door and allow the instrument to complete an inventory scan before attempting to delete the panel.

W601 Can Not Finalize

Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to finalize a panel either whose status is not Removable or that has an unignored Needs Attention flag set.	If panel is Ongoing, wait until it becomes Removable before attempting to finalize it. If panel has a Needs Attention flag set, resolve the reason before attempting to finalize the panel.

W602 Invalid Date		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a record containing, or to print a report specifying, an invalid date. Or, a panel lot barcode containing an invalid expiration date was scanned.	Enter a valid date.

W604 Panel Still Testing		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a Final ID to an ID or Combination panel that has not reported an Instrument organism ID.	A Final ID cannot be selected and saved until the instrument ID results are obtained. Wait until the instrument calculates an ID before attempting to select a Final ID.

W606 Cannot Modify Date for a Finalized Panel		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to save a panel without a panel lot number expiration date or without expiration dates for other lot number fields IF panel lot numbers have been entered.	If information for a finalized panel must be modified, recall the panel and un-finalize it by disabling the Finalize field. Make the required changes and save the information.

W607 Missing Date

Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to save a panel without a panel lot number expiration date or without expiration dates for other lot number fields IF panel lot numbers have been entered.	A valid panel lot number expiration date must be entered for panels. If other lot number fields are completed, their corresponding expiration dates must be entered.

W608 Improper Barcode Scan

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	The barcode number scanned is too long for the current active field, or is not the correct type of barcode for the field (e.g., scanning a lot number in a non-lot number field).	Verify which field is currently active and check what type of barcode you are trying to scan.

W609 Database Full

Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt was made to log in a new panel and 50 panels exist in the instrument database with a pending record status, with no panels eligible for deletion.	Pending panels should be placed in the instrument and allowed to complete an inventory scan before you attempt to log in any new panels.

W700 Invalid LIS Code		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An invalid LIS code was entered in the Organism Configuration or Antimicrobial Configuration screen.	LIS Codes can be up to 20 alphanumeric characters.

W701 Duplicate LIS Code		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	A duplicate LIS code was entered in the Organism Configuration or Antimicrobial Configuration screen.	LIS Codes must be unique in the system.

W800 Panel Lot Undefined		
Message Type	Alert Cause	Alert Corrective Action
Field Level	An attempt was made to log in a QC panel but the panel sequence number is not within the range of a defined/saved Panel Lot.	When QC Lot Support is enabled in Configuration, all QC panels must belong to a defined/ saved panel lot. Clinical panels may be from an undefined panel lot and still be logged in, though such panels generate a Panel Lot Undefined Needs Attention condition.

W801 Panel Lot Range Incomplete

Message Type	Alert Cause	Alert Corrective Action
Field Level	During Panel Lot Definition, an attempt was made to save a panel lot but one or both panel sequence number fields were blank.	Both a starting sequence number (lowest) and ending sequence number (highest) must be scanned for any panel lot to be saved. The sequence numbers can be scanned in any order, but both must be scanned. You cannot type in a sequence number on the Panel Lot Definition screen.

W802 Panel Lot Range Invalid

Message Type	Alert Cause	Alert Corrective Action
Field Level	During Panel Lot Definition, an attempt was made to save a panel lot but the two panels scanned were different panel types.	Both panels must be the same type for any panel lot to be saved.

W803 Panel Lot Range Conflict

Message Type	Alert Cause	Alert Corrective Action
Field Level	During Panel Lot Definition, an attempt was made to save a panel lot but one or both panel sequence numbers conflicted with an existing panel lot definition.	A panel sequence number can only belong to one defined/saved panel lot. Contact your local BD representative.

W903 Panel Lot Expired		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	A panel lot number barcode is scanned containing an expired expiration date.	Enter an unexpired panel log number barcode.

When Commission Duplicate Serial Number		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	There is an instrument with the same serial number already commissioned.	Enter a unique serial number.

When Commission Duplicate Instrument Identifier		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	There is an instrument with the same designation already commissioned.	Check for instrument designation duplication.

W904 Commission Invalid Serial Number		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	The serial number must be a minimum of six characters.	The serial number must be six characters.

W905 Decommission of Connected Instrument		
Message Type	Alert Cause	Alert Corrective Action
OK Print Screen	An attempt is being made to decommission an instrument.	Disconnect the instrument prior to decommissioning.

6.4 Event Log Messages

6.4.1 LIS Related Messages

The Event Log tab, located on the bottom right of the Maintenance screen enables the user to see messages that the instrument writes to the Event Log. The messages are in the following format:

date time message type: message text

where *date* represents the day, month, and year in the format you have chosen

time represents the time in the format you have chosen

message type is one of the following:

- LIS Interface Message
- LIS Unsolicited Message
- LIS Order Cancelled
- LIS Query Assembly
- LIS Receiving Query
- LIS Configuration Change

message text is the actual text message that appears

Below is a list of the messages, along with a description of the message and any actions that may be performed to correct the problem. These messages are grouped by the message type listed above.

LIS Interface Messages

LIS Interface Messages are library messages generated by the LIS manager. They are listed below in alphabetical order.

Bad Frame Received From LIS

DESCRIPTION – This error is generated when the LIS downloads a frame that has not been properly formatted. The error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – The solution to this problem is to review the information that is being sent to the instrument via a communications line monitor. Compare the information captured with the specifications found in the BD – LIS Vendor Interface Document. Any discrepancies observed should be corrected and then the transmission should be attempted again.

Detailed Download Message Attached to Notification

DESCRIPTION – These notifications are generated whenever a complete message is sent or received across the port. The error types are LIS_NOTIFY. The detailed description is not displayed.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Detailed Upload Message Attached to Notification

DESCRIPTION – These notifications are generated whenever a complete message is sent or received across the port. This provides the Host Application with the complete ASTM message string that was exchanged. The error types are LIS_NOTIFY and the detailed descriptions contain the ASTM message string.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Disallowed Characters Contained in Field

DESCRIPTION – This error is generated when an upload field contains characters that are not allowed in an upload message. These characters differ based on logical protocol and will be stripped from the upload field. This error type is LIS_NON_FATAL. The detailed description is not displayed.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Download Field Was Concatenated

DESCRIPTION – This notification is generated when a download field is larger than the maximum size set by the Host Application. The field will be concatenated to the maximum size. This error type is listed as LIS_LOG. The detailed description is not displayed.

CORRECTIVE ACTION(S) – In this case the content that is being assembled needs to be reviewed. Review the BD – Vendor Interface Specification for limitations on field lengths in the message. Then the LIS code should be updated so that the information sent to the instrument remains within the defined limits.

Download Message Has a Bad or Missing Header Record

DESCRIPTION – This error is generated when a download message does not have a properly formatted header record. The LIS IM cannot continue processing the message and it will be deleted. This error type is listed as LIS_NON_FATAL, and the detailed description holds the ASTM message received from the LIS.

CORRECTIVE ACTION(S) – Check the header line for errors. Correct the header error and resend the message.

Download Physical Communication with LIS Has Begun

DESCRIPTION – This notification is generated when the download thread in the physical layer is actively downloading a message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Download Physical Communication with LIS Has Completed

DESCRIPTION – This notification is generated when the download thread in the logical layer has completed processing a download message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Download Record is Out Of Sequence Or Has Bad Sequence Number

DESCRIPTION – This error is generated when a download record contains records with a bad sequence number, or violates the record hierarchy for the logical level protocol. This will cause the message to be deleted from the download queue. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Appropriate updates should be made to the application generating the logical message so the structure conforms to the ASTM specifications for the interface. Most questions concerning the message content can be addressed by referencing the BD – Vendor Interface Specifications and the ASTM specifications for the interface.

Early Termination of Transfer Session By LIS

DESCRIPTION – This error is generated when there is an error in LIS communication. The control character to end a transfer session was received before the appropriate number of characters for the frame were received.

CORRECTIVE ACTION(S) – Review the frames being exchanged between the instrument and the LIS to assure that the proper number of characters are contained within the packets that are being exchanged. If the correct number of characters is not in the packet, update the code so that the correct number of characters is included in each frame.

Expected Frame Not Sent

DESCRIPTION – This error is generated when there is an error in LIS communication. The LIS initiated a transfer session but did not send any data before a time-out occurred. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Connect a communications line monitor and restart the transmission. Observe the transmission to determine if the complete transmission occurs or if the transmission is interrupted early. If the transmission appears to be complete or the LIS appears to be attempting to send the transmission contact your local BD representative for assistance. If the transmission appears to be incomplete from the LIS side of the transmission then the LIS code should be investigated for problems that could terminate the transmission early.

Internal Assert Condition Found in LIS Library

DESCRIPTION – This message indicates that the Instrument has encountered an error that should not occur. This error will result in the instrument rebooting.

CORRECTIVE ACTION(S) – Contact your local BD representative. At the time the problem occurs the application writes data to the log file that indicates the nature of the Assert. The BD representative should instruct you on the appropriate procedure to collect the data for the condition.

LIS Debug Error

DESCRIPTION – This error is generated when there has been a problem internally to the instrument physical layer. This error may be triggered for an index out of bounds, or unsupported memory area IDs, etc. This error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – This error is a Fatal Error. If this error appears in the log file, contact your local BD representative. The sequence of events that produced this message should be documented so that BD representatives can reproduce the error and then provide an appropriate course of action to address the problem.

LIS Did Not Acknowledge Sent Frame

DESCRIPTION – This error is generated when the instrument has sent a frame to the LIS but has not received an acknowledgment before a time-out occurred. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – This message should not occur during normal operation of the instrument interface. If this error is being encountered during development of the interface then the LIS development group should connect a data communication monitor to the serial interface cable and review the information that is being exchanged between the two devices. It is likely that the LIS is not generating the appropriate response to the message that the instrument has sent.

LIS Is Not Responding To Output Request

DESCRIPTION – This error is generated when the instrument is trying to establish a transfer session but the LIS is not responding. When this error is sent, the instrument is assuming the LIS connection is broken or the LIS is down. This error type is listed as LIS_NON_FATAL.

CORRECTIVE ACTION(S) – Check to assure that the LIS interface is active and ready to receive messages from the instrument. If the LIS is operating correctly, connect a communications line monitor and review the information that has been captured. If the instrument is attempting to establish communications with the LIS, it would be appropriate to review the code that has been written to interface with the instrument and assure that the code is appropriate to receive the information sent by the instrument.

LIS Never Completed Current Frame

DESCRIPTION – This error is generated when there is an error in LIS communication. The control characters expected to end a transmitted frame were never received.

CORRECTIVE ACTION(S) – Review the frames being exchanged between the instrument and the LIS via line monitor. Be sure to review all frames included in the transmission. If any frames do not contain the appropriate termination characters make appropriate changes to the interface code to correct the problem.

Logical Processing of LIS Data Has Begun

DESCRIPTION/CORRECTIVE ACTION(S) – Message is informational – no action required.

Logical Processing of LIS Data Has Completed

DESCRIPTION/CORRECTIVE ACTION(S) – Message is informational – no action required.

Message Packet Passed to Host Application

DESCRIPTION – This notification is generated when the LIS IM has successfully passed a message on to the Host Application. This notification exists mostly as a debug message. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – No action required by user. Information is included to help LIS manufacturers debug and implement their interface.

Message Received By LIS

DESCRIPTION – This notification is generated after the LIS has received and acknowledged a message from the instrument.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Message Received From LIS

DESCRIPTION – This notification is generated when a download message has been properly received by the Physical Interface from the LIS. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – This message is an informational message and requires no action to be taken. The transmission was successful. This message will be generated even if the message contained content errors.

Must Re-send Output Frame to LIS

DESCRIPTION – This notification is generated when a frame sent by the LIS IM was not properly received by the LIS. The message frame will be resent, according to low level protocol specifications. This error type is listed as LIS_LOG.

CORRECTIVE ACTION(S) – If this message appears infrequently and the messages that are being exchanged are completed correctly, it is likely that no action is required. If the message above is encountered frequently, there is likely a problem in the interface. The LIS manufacturer and BD representative should be contacted to diagnose the interface to assure that it is operating properly. When the link is operating optimally this message should not appear in the log.

No Queue Memory for Download Messages

DESCRIPTION – The BD Phoenix M50 Instrument is designed to operate in a limited amount of memory (as defined in the configuration structure). When these memory resources are full with download or upload messages, this notification will be generated. The notification may be common if the LIS tries to download too much information at once. The BD Phoenix M50 Instrument will NAK data for which it cannot allocate memory, and the LIS will have to resend the data.

CORRECTIVE ACTION(S) – Resend the messages that were rejected by the instrument.

No Queue Memory for Upload Messages

DESCRIPTION – The BD Phoenix M50 Instrument is designed to operate in a limited amount of memory (as defined in the configuration structure). When these memory resources are full with download or upload messages, this notification will be generated. This message should not be encountered when using the instrument LIS interface.

CORRECTIVE ACTION(S) – If this message is encountered, the problem should be reported to the BD representative. It will be helpful if the process that generated the problem is documented.

No Response Received From Previous LIS Request Message

DESCRIPTION – This notification is generated when the LIS IM has uploaded a query but no response is received. This notification indicates that the original query was cancelled. This notification should occur some time after the query, as determined by the value in the configuration structure. The error type is listed as LIS_LOG. The detailed description is not displayed in the current version.

CORRECTIVE ACTION(S) – Validate that the LIS system is correctly connected to the instrument and that the LIS interface has been activated.

Operating System Error

DESCRIPTION – This error is generated when the operating system class encounters an error in one of its routines. This could be caused by a number of OS errors, including not properly initializing a port or not properly creating an event handle, and others. These errors should not occur under normal conditions and the LIS IM cannot recover from them. The error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Output Message Was Sent To LIS

DESCRIPTION – This notification is generated when an upload message has been successfully transmitted to the LIS. This error type is listed as LIS_NOTIFY. Notification sent when the instrument begins sending a message to the LIS. Receiving this notification does not signify the message was accepted by the LIS.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Queue Memory for Download Messages Free

DESCRIPTION – This notification is used in pairs with the No Queue Memory For Download Messages above. When memory has previously been determined to be full, and now has been released, this notification will be generated.

CORRECTIVE ACTION(S) – This is an indication to a user of a batch-oriented interface that the instrument is capable of receiving the next group of messages.

Queue Memory for Upload Messages Free

DESCRIPTION – This notification is used in pairs with the No Queue Memory For Upload Messages above. When memory has previously been determined to be full, and now has been released, this notification will be generated.

CORRECTIVE ACTION(S) – This message should not appear in the log during operation of the interface. If this message appears in the log, contact your local BD representative.

Response Message Received

DESCRIPTION – This notification is generated when the LIS has downloaded a response to a BD generated query. This is also an indication that the LIS IM can upload another query to the LIS. The error type is LIS_NOTIFY. The detailed description is not displayed in the current version.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Unsupported Field in Configuration File

DESCRIPTION – This error is generated during startup when a field in the configuration structure does not match one of the supported configurations of the LIS IM. This error type is listed as LIS_FATAL.

CORRECTIVE ACTION(S) – If this error is reported the problem should be reported to your local BD representative for further investigation.

Upload Physical Communication With LIS Has Begun

DESCRIPTION – This notification is generated when the upload thread in the physical layer is actively uploading a message.

CORRECTIVE ACTION(S) – Message is informational – no action required.

Upload Physical Communication With LIS Has Completed

DESCRIPTION – This notification is generated when the upload thread in the physical layer has completed sending a message to the LIS.

CORRECTIVE ACTION(S) – Message is informational – no action required.

LIS Unsolicited Message

The following sections contain the possible messages associated with Event Log entries in response to unsolicited download requests. They are listed below in alphabetical order.

LIS Response Error

DESCRIPTION – This message is set after the LIS sends an incorrectly formatted query message to the instrument. We will respond by sending a response error message to the LIS.

CORRECTIVE ACTION(S) – If the interface is operating correctly, this message should not be present in the log. If this message is found in the log during validation or development of the interface, the code responsible for generating the queries should be reviewed and adjusted so that the query is assembled without the invalid field. The invalid field will be the last "LIS Query Assembly" message listed before this message in the log file.

LIS Sent Query Before Previous Query Completed

DESCRIPTION – This message is set if the LIS initiates a query and while the instrument is busy with this query, the LIS cancels the first query and immediately starts another.

CORRECTIVE ACTION(S) – If the interface is operating correctly this message should not be present in the log. If this message is found in the log during validation or development of the interface, the code responsible for generating the queries should be reviewed and adjusted so that it will not send another query while it is waiting for an outstanding query to be canceled.

LIS Order Canceled Messages

DESCRIPTION –

•

The validation for the noted field failed.

Invalid Sequence Number Field

Missing Accession Number Field

Invalid Accession Number Field [cannot be 12 digits and begin with 42 or 50–59; see Section 3.3 Panel Login Tab for other requirements.

Missing Isolate Field

Invalid Isolate Field

Invalid Panel Usage

Invalid Test ID Field

Invalid Test Strain for this QC Panel Type

Mismatch Sequence Number and Test ID

Missing Test ID Field

Invalid Organism ID Field

Invalid Priority Field

Problem Storing Record to Database

When an error is detected, the message for that error is logged and the remainder of the message checking is terminated. There could be other errors in the order that were not reported. The error checking priority depends upon whether the order contains a sequence number or not.

CORRECTIVE ACTION(S) – In an operational interface it is possible to encounter these errors if the LIS implementation allows entry of data into the fields that violate the instrument field rules. In this case the user of the LIS should limit the characters entered into the fields sent to the instrument to those characters noted in the System User's Manual as valid for the field. Once the information has been updated to conform with the instrument's rules the order should be sent to the instrument again.

If the interface is under development the developer of the interface should limit the data being sent to the instrument to those characters that are appropriate for the fields.

LIS Query Assembly Messages

All of the following messages are written to the log file regardless of whether they are invalid or valid. They are placed in the log so that they can be used in conjunction with the messages that follow them. If an error is found in one of the fields, the contents of the field in error will be logged.

(Message consists of a field name shown below at **.)

DESCRIPTION –	The message is logged when the instrument uploads a query to the LIS.
**	Accession Number
	Sequence Number
	Test ID
	Test Status
	Result Qualifier
	Time Qualifier
	Start Date/Time
	End Date/Time

CORRECTIVE ACTION(S) – If the query is valid, the message is informational – no action required. If the query is invalid, then review the field in error, correct, and resend the query.

(Message consists of a field name and value shown below at ***)

DESCRIPTION – This message is logged when the instrument uploads a query to the LIS.
*** Sequence Number and the value

CORRECTIVE ACTION(S) – Message is informational – no action required.

LIS Configuration Change Messages

(Message consists of a configuration value shown below at ** with old and new values.)

DESCRIPTION –	The message is logged when the LIS configuration is changed.
**	LIS Enabled Value Changed
	Send Interpretation Results Option Changed
	Unsolicited Queries Option Changed
	Send When Placed In Instrument Option Changed
	Results Upload Options Changed
	Baud Value Changed
	Data Bits Value Changed
	Parity Value Changed
	Stop Bits Value Changed
	Packed Frames Value Changed

CORRECTIVE ACTION(S) – Message is informational – no action required.

7 – Performance, Limitations, Organism Codes, and Panel Information

7.1 Performance Characteristics

Definitions:

Essential Agreement (EA): Essential agreement occurs when the MIC of the BD Phoenix system and the CLSI Reference Broth Microdilution are identical or within ± 1 dilution of each other.

Category Agreement (CA): Category agreement occurs when the BD Phoenix system results are within ± 1 dilution from the CLSI Reference Broth Microdilution with respect to the CLSI categorical interpretative criteria.

7.1.1 Gram Negative Performance

Gram Negative Identification

In two internal studies, the performance of the BD Phoenix Gram Negative Identification was evaluated. The 0.5 inoculum density configuration and the 0.25 inoculum density configuration were tested with 721 strains (0.5) and 784 strains (0.25) respectively. Enteric and non-enteric isolate results were evaluated against commercial and noncommercial methods.

The BD Phoenix Gram Negative Identification performance is outlined below:

	Inoculum Density (McFarland)	Agreement	No Agreement	No ID
Genus Level	0.5	97.3%	2.1%	0.6%
	0.25	98.5%	1.0%	0.5%
Species Level	0.5	95.6%	3.6%	0.8%
	0.25	98.1%	1.4%	0.5%

Gram Negative Susceptibility

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Amikacin	AN	974	94.3	974	99.2
Amoxicillin	AMX	725	95.9	839	96.9
Amoxicillin-Clavulanate	AMC	636	97.8	749	97.1
Amoxicillin-Clavulanate (f)	AXC	364	96.4	411	97.6
Ampicillin	AM	639	97.5	753	98.0
Ampicillin-Sulbactam	SAM	848	96.8	962	96.9
Ampicillin-Sulbactam (f)	SXA	1682	92.3	1682	94.5
Arbekacin	ARB	2083	94.1	2083	99.2
Aztreonam	ATM	1488	97.8	1431	98.5
Cefazolin	CZ	634	97.0	752	97.1
Cefdinir	CDR	633	92.7	747	94.5
Cefditoren	CDN	1270	97.2	1270	98.4
Cefepime	FEP	1463	97.8	1463	98.6
Cefetamet	CAT	629	96.0	746	97.7
Cefixime	CFM	1454	96.7	1454	98.1
Cefmetazole	CMZ	608	94.7	718	96.4
Cefoperazone	CFP	854	95.1	972	97.7
Cefoperazone-Sulbactam	SCP	1921	93.1	1921	96.0
Cefotaxime	CTX	849	96.9	970	97.2
Cefotetan	CTT	627	96.2	748	96.4
Cefotiam	CFT	1428	94.9	1428	97.6
Cefoxitin	FOX	628	96.7	748	97.6
Cefpirome	CPO	846	96.2	964	97.0
Cefpodoxime	CPD	616	94.5	737	97.0
Cefsulodin	CFS	55	98.2	56	98.2
Ceftazidime	CAZ	848	95.9	969	97.5
Ceftibuten	CTB	590	90.5	708	95.9
Ceftizoxime	ZOX	854	97.3	971	97.2
Ceftriaxone	CRO	2416	96.1	2416	91.6
Cefuroxime	CXM	623	97.0	744	98.5
Cephalexin	CN	319	97.8	366	96.5

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Cephalexin	CF	613	98.0	731	98.5
Chloramphenicol	C	978	97.7	978	98.7
Ciprofloxacin	CIP	977	99.4	977	99.7
Colistin	CL	467	97.9	467	98.7
Ertapenem*	ETP	1632	97.7	1265	99.4
Fosfomycin	FF	446	93.0	446	96.9
Garenoxacin	GRN	1977	98.1	1977	99.4
Gatifloxacin	GAT	752	99.5	752	99.6
Gemifloxacin	GEM	2096	98.6	2096	99.0
Gentamicin	GM	973	95.8	973	99.8
Imipenem	IPM	1545	94.8	1483	99.0
Isepamicin	ISP	468	93.8	468	99.4
Kanamycin	K	735	95.8	735	99.9
Levofloxacin	LVX	972	99.6	972	100.0
Lomefloxacin	LOM	976	99.1	976	99.6
Mecillinam	MEC	345	95.7	345	94.8
Meropenem	MEM	1450	97.0	1394	98.3
Minocycline	MI	2094	93.8	2094	98.3
Moxalactam	MOX	2063	96.9	2063	98.3
Moxifloxacin	MXF	746	98.5	747	99.7
Nalidixic Acid	NA	750	94.0	750	99.1
Netilmicin	NET	974	96.7	974	99.3
Nitrofurantoin	FM	744	98.4	744	98.8
Norfloxacin	NOR	976	98.9	976	99.6
Ofloxacin	OFX	971	99.3	971	99.7
Pefloxacin	PEF	469	98.5	469	99.4
Piperacillin	PIP	860	94.7	973	97.1
Piperacillin-Tazobactam	TZP	856	92.9	970	96.2
Temocillin	TEM	1410	96.1	1410	98.9
Tetracycline	TE	975	95.6	975	98.5
Ticarcillin	TIC	859	94.8	973	97.8
Ticarcillin-Clavulanate	TIM	534	92.7	589	94.6
Tigecycline	TGC	1428	97.1	1110	96.7
Tobramycin	NN	977	94.2	977	99.4
Trimethoprim	TMP	752	96.0	752	99.6
Trimeth-Sulfa (DIN)	STG	463	97.8	463	97.8
Trimethoprim-Sulfamethoxazole	SXT	976	95.5	976	97.3

* Ability to Detect Resistance Unknown

7.1.2 Gram Positive Performance

Gram Positive Identification

In two internal studies, the performance of the BD Phoenix Gram Positive Identification was evaluated. The 0.5 inoculum density configuration and the 0.25 inoculum density configuration were tested with 696 strains (0.5) and 755 strains (0.25) respectively. Results were evaluated against commercial and non-commercial methods.

The BD Phoenix Gram Positive Identification performance is outlined below:

	Inoculum Density (McFarland)	Agreement	No Agreement	No ID
Genus Level	0.5	98.1%	0.3%	0.6%
	0.25	99.6%	0.0%	0.4%
Species Level	0.5	95.4%	3.9%	0.7%
	0.25	98.0%	1.6%	0.4%

Gram Positive Susceptibility

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Amikacin	AN	487	95.7	487	95.7
Amoxicillin	AMX	395	90.6	659	97.1
Amoxicillin-Clavulanate	AMC	180	97.8	446	97.5
Amoxicillin-Clavulanate (f)	AXC	397	96.7	664	97.6
Ampicillin	AM	402	94.0	667	98.7
Ampicillin-Sulbactam	SAM	179	96.7	445	97.1
Ampicillin-Sulbactam (f)	SXA	1449	94.6	1449	95.2
Arbekacin	ARB	973	96.2	973	99.9
Azithromycin	AZM	702	95.7	702	97.7
Cefaclor	CEC	182	91.8	449	95.6
Cefazolin	CZ	180	97.2	441	97.5
Cefdinir	CDR	182	94.5	434	96.8
Cefditoren	CDN	944	96.2	944	97.8
Cefepime	FEP	181	97.2	446	97.5
Cefmetazole	CMZ	174	96.6	430	97.2
Cefoperazone	CFP	184	95.7	447	97.8
Cefotaxime	CTX	185	96.2	446	97.8
Cefotetan	CTT	183	94.5	439	96.6
Cefotiam	CFT	965	91.5	965	94.0
Cefoxitin	FOX	184	96.2	445	97.5
Cefozopran	CFZ	1130	89.8	1130	89.6
Cefpirome	CPO	183	100.0	446	97.5
Cefpodoxime	CPD	185	91.4	451	95.6

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Ceftaroline*	CPT	1313	94.5	866	99.7
Ceftazidime	CAZ	183	95.1	448	97.1
Ceftizoxime	ZOX	184	93.5	447	96.2
Ceftriaxone	CRO	184	95.1	447	97.3
Cefuroxime	CXM	184	94.6	447	98.2
Cephalexin	CN	185	90.3	452	97.4
Cephalothin	CF	182	97.8	441	97.3
Chloramphenicol	C	705	92.8	705	95.3
Ciprofloxacin	CIP	686	96.4	686	97.4
Clarithromycin	CLR	485	96.3	485	98.1
Clindamycin	CC	703	98.2	703	98.3
Daptomycin*	DAP	1361	97.6	1361	99.8
Doxycycline	D	1211	96.5	1211	99.8
Ertapenem	ETP	670	91.6	670	95.4
Erythromycin	E	474	95.4	474	95.8
Fosfomycin	FF	472	97.0	472	98.3
Fusidic Acid	FA	477	96.7	477	99.0
Garenoxacin	GRN	1212	97.1	1212	99.0
Gatifloxacin	GAT	481	98.8	481	100.0
Gemifloxacin	GEM	1400	96.2	1400	97.4
Gentamicin	GM	487	98.4	487	99.2
Gentamicin-Syn	GMS	NA	NA	198	97.5
Gentamicin-Syn (SFM)	GMF	NA	NA	87	100.0
Imipenem	IPM	185	100.0	448	97.5
Inducible Macrolide Resistance (iMLSB) Test	ECC	NA	NA	295	97.6
Kanamycin	K	463	92.2	464	97.4
Kanamycin Synergy	KS	501	100.0	501	96.6
Levofloxacin	LVX	698	96.9	698	98.6
Lincomycin	L	696	98.3	696	98.7
Linezolid	LZD	474	99.4	474	100.0
Lomefloxacin	LOM	486	97.7	486	99.2
Meropenem	MEM	185	98.4	447	97.5
Minocycline	MI	1448	98.5	1448	99.2
Moxalactam	MOX	935	95.9	935	96.0
Moxifloxacin	MXF	486	99.0	486	99.6
Mupirocin	MUP	938	96.6	938	99.7
Mupirocin High Level	MUH	968	100.0	968	99.9

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Netilmicin	NET	488	93.9	488	96.3
Nitrofurantoin	FM	707	98.3	707	98.3
Norfloxacin	NOR	674	96.0	674	97.8
Ofloxacin	OFX	487	99.2	488	99.8
Oxacillin	OX	449	95.8	449	97.8
Pefloxacin	PEF	706	97.9	706	98.0
Penicillin	P	401	96.0	662	98.0
Piperacillin	PIP	382	81.7	639	98.0
Piperacillin-Tazobactam	TZP	386	100.0	649	98.3
Pristinamycin	PR	695	93.2	695	95.5
Quinupristin-Dalfopristin	SYN	704	97.2	704	97.6
Rifampin	RA	700	84.0	700	88.3
Streptomycin-Syn	STS	NA	NA	198	99.0
Streptomycin-Syn (SFM)	STF	NA	NA	87	98.8
Teicoplanin	TEC	671	97.6	671	99.4
Telithromycin	TEL	898	97.3	898	98.9
Tetracycline	TE	707	96.3	707	98.6
Ticarcillin	TIC	395	95.4	659	98.2
Ticarcillin-Clavulanate	TIM	396	97.5	656	97.9
Tigecycline*	TGC	1356	96.8	1044	100
Tobramycin	NN	473	97.0	473	97.3
Trimethoprim	TMP	479	95.4	479	97.9
Trimeth-Sulfa (DIN)	STG	645	91.0	646	92.7
Trimethoprim-Sulfamethoxazole	SXT	654	89.5	654	93.0
Vancomycin	VA	702	99.3	702	99.9

* Ability to Detect Resistance Unknown

7.1.3 Streptococci Performance (with BD Phoenix SMIC/ID, SMIC Panels)

Streptococci Identification (SMIC/ID only)

In an internal study, the performance of the BD Phoenix Streptococci Identification was evaluated. Results from 655 isolates were evaluated against commercial and noncommercial methods.

The BD Phoenix Streptococci Identification performance is outlined below:

	Agreement	No Agreement	No ID
Species Level	96.3%	2.4%	1.2%

Streptococci Susceptibility

The performance of the BD Phoenix Streptococci AST System was evaluated at four clinical trial sites using clinical isolates. Comparisons were made to AST results generated from reference broth microdilution panels prepared according to CLSI standard guidelines. Discrepant results were arbitrated by duplicate repeat testing in both systems.

ANTIMICROBIAL	CODE	EA N	EA %	CA N	CA %
Amoxicillin	AMX	1576	99.2	1576	99.7
Amoxicillin-Clavulanate	AMC	1564	97.9	1564	99.3
Ampicillin	AM	1569	98.0	1569	99.4
Cefepime	FEP	1571	97.6	1571	99.7
Cefotaxime	CTX	1578	98.5	1578	99.6
Ceftriaxone	CRO	1579	98.5	1579	99.8
Cefuroxime	CXM	1581	97.7	1581	98.8
Chloramphenicol	C	1587	97.4	1587	99.6
Clindamycin	CC	1591	94.0	1591	97.4
Daptomycin	DAP	1566	93.1	1566	99.9
Ertapenem*	ETP	1585	97.9	1585	99.6
Erythromycin	E	1577	93.2	1577	98.2
Garenoxacin*	GRN	1515	95.5	1515	99.7
Gatifloxacin	GAT	1587	95.1	1587	99.5
Gemifloxacin*	GEM	1592	98.9	1592	99.7
Gentamicin-Syn	GMS	1032	99.6	1032	99.6
Imipenem	IPM	1581	98.0	1581	99.6
Kanamycin-Syn	KS	1032	99.6	1032	99.6
Levofloxacin	LVX	1595	97.4	1595	99.8
Linezolid	LZD	1586	97.1	1586	99.4
Meropenem	MEM	1579	97.4	1579	99.7
Moxifloxacin*	MXF	1590	95.8	1590	99.6
Ofloxacin	OFX	1594	98.2	1594	98.7
Penicillin	P	1585	93.9	1585	98.5
Pristinamycin*	PR	1583	94.8	1583	99.8
Quinupristin-dalfopristin*	SYN	1581	97.3	1581	99.7
Streptomycin-Syn	STS	1031	99.9	1031	99.9
Teicoplanin*	TEC	1593	99.8	1593	99.9
Telithromycin*	TEL	1574	95.9	1574	98.5
Tetracycline	TE	1589	93.3	1589	97.6
Trimethoprim-Sulfamethoxazole	SXT	1008	93.8	1008	98.3
Vancomycin*	VA	1588	98.2	1588	99.9

* Ability to Detect Resistance Unknown

7.1.4 Identification of Yeast Species

The performance of the BD Phoenix Yeast identification was evaluated across multiple sites using pure colonies isolated from Sabouraud Dextrose Agar (SAB) and BD Trypticase Soy Agar with 5% Sheep Blood (TSA). Results from 519 (SAB) and 510 (TSA) clinical and challenge isolates were evaluated against conventional and molecular methods.

The BD Phoenix Yeast identification performance is outlined below:

	Source Media	Agreement	No Agreement	No ID
Genus/Species Level	SAB	95.2%	3.8%	1.0%
	TSA	96.5%	2.7%	0.8%

Additionally, testing was performed at multiple sites to demonstrate reproducibility. The identification results obtained using the BD Phoenix system were compared with expected results. This performance testing demonstrated inter-site reproducibility of $\geq 95\%$.

7.2 Limitations of the Procedure

See the package insert shipped with the panel for specific organism/antimicrobial limitations.

General

A Gram stain test is required for the selection of the appropriate BD Phoenix panel types. Accurate identification and/or AST results may not be made without this test.

Use only well-isolated bacterial or yeast colonies from one of the recommended primary isolation media. Use of mixed colonies could result in inaccurate identification and/or AST interpretations.

If the instrument inoculum density (for the panel type being used) is configured to 0.5, an inoculum density of 0.50–0.60 McFarland must be met. Only the BD BBL CrystalSpec or BD PhoenixSpec Nephelometer can be used to measure the inoculum density.

If the instrument inoculum density (for the panel type being used) is configured to 0.25, an inoculum density of 0.20–0.30 McFarland must be met. Only the BD PhoenixSpec Nephelometer can be used to measure the inoculum density for this range.

For identification of yeast, a suspension equivalent of 2.00–2.40 McFarland standard must be met and prepared only with the BD PhoenixSpec Nephelometer. Use of alternate methods for suspension preparation may cause erroneous identification results.

BD Phoenix panels can be read only by the BD Phoenix M50 instrument. Visual interpretation of the BD Phoenix panels is not possible. Any attempt to manually interpret results from the panel may lead to misidentification and/or inaccurate AST interpretations.

Panels must be placed into the BD Phoenix M50 instrument within 30 minutes of inoculation.

For the most reliable results, it is recommended that the QC organisms be subcultured at least twice on two consecutive days onto TSA with 5% Sheep Blood before use in the BD Phoenix system.

Identification

The unique panel environment combined with the shortened incubation time may result in BD Phoenix panel reactions varying from those obtained using conventional biochemical media.

Antimicrobial Susceptibility Testing

AST Indicator solution may be added to the AST broth tubes and stored in the dark for up to 8 hours prior to use. Use only one free-falling drop of the AST Indicator solution from the dropper bottle. If more than one drop is added inadvertently, discard the tube and use a fresh tube of AST broth.

After the addition of BD Phoenix AST Indicator Solution to the AST broth tubes, mix by inversion. DO NOT VORTEX. Vortexing may cause air bubbles to form in the AST broth, which can result in inappropriate filling of the BD Phoenix panel during inoculation.

Because of the low probability of occurrence some organisms included in the ID taxa are not included in the AST database. These organisms will display the message Organism not included in the AST database, perform alternate method.

For some organism/antimicrobial combinations, the absence of resistant strains precludes defining any result categories other than susceptible. For strains yielding results suggestive of a nonsusceptible category, organism identification and antimicrobial susceptibility test results should be confirmed. Subsequently, the isolates should be saved and submitted to a reference laboratory that will confirm the result using the CLSI reference microbroth dilution method.

The use of CHROMagar Orientation may produce false susceptibility results when testing erythromycin with gram-positive organisms. Antimicrobial susceptibility results should be confirmed using BD Trypticase Soy Agar with 5% Sheep Blood.

This media type [Chocolate Agar] should not be used for Streptococcal identification with SMIC/ID panels. Chocolate Agar may be used for Streptococcal susceptibility testing only.

The use of Columbia Agar with 5% Horse Blood may produce significantly higher MIC for SXT with *Streptococcus* species, which may result in false resistance. Antimicrobial susceptibility test results should be confirmed using BD Trypticase Soy Agar with 5% Sheep Blood.

NOTES

MIC results for isolates without a genus and species identified in the BD Phoenix results Final ID field may not be valid.

When changing the BD Phoenix Final ID from an organism not claimed in the AST Taxa list to an organism that is claimed in the AST Taxa, the MIC and Interpretation results will be based on the Final ID that is claimed in the AST Taxa.

7.3 Organism Codes, Short Names, Long Names

Long Name	Short Name	Code
<i>Achromobacter denitrificans</i>	Achr. denitrificans	ALCDEN
<i>Achromobacter piechaudii</i>	Achr. piechaudii	ALCPIE
<i>Achromobacter species</i>	Achr. species	ACHRSPE
<i>Achromobacter xylosoxidans</i>	Achr. xylosoxidans	ALCXYL

Long Name	Short Name	Code
<i>Acinetobacter baumannii</i>	Acinet. baumannii	ACINBAU
<i>Acinetobacter baumannii/calcoaceticus</i> complex	Acinet. baumannii/calco. cplx	ACINBCX
<i>Acinetobacter baumannii/haemolyticus</i>	Acinet. baumannii/haemolyticus	ACINBAUHAE
<i>Acinetobacter calcoaceticus</i>	Acinet. calcoaceticus	ACINCAL
<i>Acinetobacter haemolyticus</i>	Acinet. haemolyticus	ACINHAE
<i>Acinetobacter johnsonii</i>	Acinet. johnsonii	ACINJOH
<i>Acinetobacter junii</i>	Acinet. junii	ACINJUN
<i>Acinetobacter lwoffii</i>	Acinet. lwoffii	ACINLWO
<i>Acinetobacter lwoffii/haemolyticus</i>	Acinet. lwoffii/haemol.	ACINLWOHAE
<i>Acinetobacter radioresistens</i>	Acinet. radioresistens	ACINRAD
<i>Acinetobacter species</i>	Acinet. species	ACINSPE
<i>Actinobacillus lignieresii</i>	Actinob. lignieresii	ACTBLIG
<i>Actinobacillus suis</i>	Actinob. suis	ACTBSUI
<i>Actinobacillus ureae</i>	Actinob. ureae	ACTBURE
<i>Aerococcus species</i>	Aeroc. species	AERCSPE
<i>Aerococcus urinae</i>	Aeroc. urinae	AERCURI
<i>Aerococcus viridans</i>	Aeroc. viridans	AERCVIR
<i>Aeromonas allosaccharophila</i>	Aerom. allosaccharophila	AERMALL
<i>Aeromonas caviae</i>	Aerom. caviae	AERMCAV
<i>Aeromonas eucrenophila</i>	Aerom. eucrenophila	AERMEUC
<i>Aeromonas hydrophila</i>	Aerom. hydrophila	AERMHYD
<i>Aeromonas hydrophila</i> group	Aerom. hydrophila gr.	AERMHYDGR
<i>Aeromonas jandaei</i>	Aerom. jandaei	AERMJAN
<i>Aeromonas media</i>	Aerom. media	AERMMED
<i>Aeromonas salmonicida</i>	Aerom. salmonicida	AERMSAL

Long Name	Short Name	Code
<i>Aeromonas salmonicida</i> ssp <i>achromogenes</i>	Aerom. salmonic. ssp ach.	AERMSALA
<i>Aeromonas salmonicida</i> ssp <i>masoucida</i>	Aerom. salmonic. ssp mas.	AERMSALM
<i>Aeromonas salmonicida</i> ssp <i>pectinolytica</i>	Aerom. salmonic. ssp pec.	AERMSALPE
<i>Aeromonas salmonicida</i> ssp <i>salmonicida</i>	Aerom. salmonic. ssp sal.	AERMSALSA
<i>Aeromonas salmonicida</i> ssp <i>smithia</i>	Aerom. salmonic. ssp smit.	AERMSALSM
<i>Aeromonas schubertii</i>	Aerom. schubertii	AERMSCH
<i>Aeromonas</i> species	Aerom. species	AERMSPE
<i>Aeromonas trota</i>	Aerom. trota	AERMTRO
<i>Aeromonas veronii</i> bv <i>sobria</i>	Aerom. veronii bv sobria	AERMVERS
<i>Aeromonas veronii</i> bv <i>veronii</i>	Aerom. veronii bv veronii	AERMVERV
<i>Alcaligenes faecalis</i>	Alc. faecalis	ALCFAE
<i>Alcaligenes faecalis</i> ssp <i>faecalis</i>	Alc. faecalis ssp faecalis	ALCFAEF
<i>Alcaligenes</i> species	Alc. species	ALCSPE
<i>Alloiococcus otitis</i>	All. otitis	ALLOTI
<i>Arcanobacterium haemolyticum</i>	Arcan. haemolyticum	ARCAHAE
<i>Bacillus cereus</i>	Baci. cereus	BACICER
<i>Bacillus circulans</i>	Baci. circulans	BACICIR
<i>Bacillus coagulans</i>	Baci. coagulans	BACICOA
<i>Bacillus licheniformis</i>	Baci. licheniformis	BACILIC
<i>Bacillus megaterium</i>	Baci. megaterium	BACIMEG
<i>Bacillus pumilus</i>	Baci. pumilus	BACIPUM
<i>Bacillus subtilis</i>	Baci. subtilis	BACISUB
<i>Bacillus thuringiensis</i>	Baci. thuringiensis	BACITHU
<i>Bergeyella zoohelcum</i>	Ber. zoohelcum	BERZOO
<i>Bordetella bronchiseptica</i>	Bord. bronchiseptica	BORBROS
<i>Brevibacillus brevis</i>	Brevs. brevis	BACIBRE

Long Name	Short Name	Code
<i>Brevibacterium</i> species	Brevm. species	BREISPE
<i>Brevundimonas diminuta</i>	Brevu. diminuta	BREUDIM
<i>Brevundimonas</i> species	Brevu. species	BREUSPE
<i>Brevundimonas vesicularis</i>	Brevu. vesicularis	BREUVES
<i>Burkholderia caryophylli</i>	Burk. caryophylli	BURCAR
<i>Burkholderia cepacia</i> complex	Burk. cepacia complex	BURCEP
<i>Burkholderia cepacia</i> / <i>Ralstonia pickettii</i>	Burk. cepacia/Ral. pickettii	BURCEPRALPIC
<i>Burkholderia gladioli</i>	Burk. gladioli	BURGLA
<i>Burkholderia glathei</i>	Burk. glathei	BURGLT
<i>Burkholderia graminis</i>	Burk. graminis	BURGRA
<i>Burkholderia multivorans</i>	Burk. multivorans	BURMUL
<i>Burkholderia phenazinium</i>	Burk. phenazinium	BURPHE
<i>Burkholderia pyrrocinia</i>	Burk. pyrrocinia	BURPYR
<i>Burkholderia</i> species	Burk. species	BURSPE
<i>Burkholderia</i> species/ <i>Ralstonia</i> species	Burk. species/Ral. species	BURSPERALSPE
<i>Candida albicans</i>	Can. albicans	CANALB
<i>Candida apicola</i>	Can. apicola	CANAPI
<i>Candida boidinii</i>	Can. boidinii	CANBOI
<i>Candida bracarensis</i>	Can. bracarensis	CANBRA
<i>Candida catenulata</i>	Can. catenulata	CANCAT
<i>Candida ciferrii</i>	Can. ciferrii	CANCIF
<i>Candida dubliniensis</i>	Can. dubliniensis	CANDUB
<i>Candida firmetaria</i>	Can. firmetaria	CANLAM
<i>Candida freyschussii</i>	Can. freyschussii	CANFRE
<i>Candida glabrata</i>	Can. glabrata	TORGLA
<i>Candida guilliermondii</i>	Can. guilliermondii	CANGUI

Long Name	Short Name	Code
<i>Candida guilliermondii</i> var <i>membranaefaciens</i>	Can. guillier. var membranaef.	CANGUIM
<i>Candida haemulonii</i>	Can. haemulonii	CANHAE
<i>Candida inconspicua</i>	Can. inconspicua	CANINC
<i>Candida kefyr</i>	Can. kefyr	CANKEF
<i>Candida krusei</i>	Can. krusei	CANKRU
<i>Candida lipolytica</i>	Can. lipolytica	CANLIP
<i>Candida lusitanae</i>	Can. lusitanae	CANLUS
<i>Candida magnoliae</i>	Can. magnoliae	CANMAG
<i>Candida melibiosica</i>	Can. melibiosica	CANMEL
<i>Candida membranifaciens</i>	Can. membranifaciens	CANMEM
<i>Candida norvegensis</i>	Can. norvegensis	CANNOR
<i>Candida parapsilosis</i> complex	Can. parapsilosis complex	CANPARPX
<i>Candida pararugosa</i>	Can. pararugosa	CANPARR
<i>Candida pelliculosa</i>	Can. pelliculosa	CANPEL
<i>Candida pulcherrima</i>	Can. pulcherrima	CANPUL
<i>Candida rugosa</i>	Can. rugosa	CANRUG
<i>Candida sake</i>	Can. sake	CANSAK
<i>Candida sphaerica</i>	Can. sphaerica	CANSPH
<i>Candida tropicalis</i>	Can. tropicalis	CANTRO
<i>Candida utilis</i>	Can. utilis	CANUTI
<i>Candida viswanathii</i>	Can. viswanathii	CANVIS
<i>Candida zeylanoides</i>	Can. zeylanoides	CANZEY
<i>Cardiobacterium hominis</i>	Card. hominis	CARHOM
CDC group Vb-3	CDC Vb-3	CDCVb3
<i>Cedecea davisae</i>	Ced. davisae	CEDDAV
<i>Cedecea lapagei</i>	Ced. lapagei	CEDLAP

Long Name	Short Name	Code
<i>Cedecea neteri</i>	Ced. neteri	CEDNET
<i>Cedecea species</i>	Ced. species	CEDSPE
<i>Cedecea species 3</i>	Ced. species 3	CEDSPE3
<i>Cedecea species 5</i>	Ced. species 5	CEDSPE5
<i>Cellulomonas turbata</i>	Cell. turbata	OERTUR
<i>Cellulosimicrobium cellulans</i>	Cellulo. cellulans	OERXAN
<i>Chromobacterium violaceum</i>	Chrom. violaceum	CHROVIO
<i>Chryseobacterium gleum</i>	Chryseob. gleum	CHRBGLE
<i>Chryseobacterium indologenes</i>	Chryseob. indologenes	CHRBIND
<i>Chryseobacterium scophthalmum</i>	Chryseob. scophthalmum	CHRBSCO
<i>Chryseobacterium species</i>	Chryseob. species	CHRBSP
<i>Citrobacter amalonaticus</i>	Cit. amalonaticus	CITAMA
<i>Citrobacter braakii</i>	Cit. braakii	CITBRA
<i>Citrobacter farmeri</i>	Cit. farmeri	CITFAR
<i>Citrobacter freundii</i>	Cit. freundii	CITFRE
<i>Citrobacter gillenii</i>	Cit. gillenii	CITSPE10
<i>Citrobacter koseri</i>	Cit. koseri	CITKOS
<i>Citrobacter murlinae</i>	Cit. murlinae	CITSPE11
<i>Citrobacter rodentium</i>	Cit. rodentium	CITSPE9
<i>Citrobacter sedlakii</i>	Cit. sedlakii	CITSED
<i>Citrobacter species</i>	Cit. species	CITSPE
<i>Citrobacter werkmanii</i>	Cit. werkmanii	CITWER
<i>Citrobacter youngae</i>	Cit. youngae	CITYOU
<i>Comamonas terrigena</i>	Coma. terrigena	COMTER
<i>Comamonas testosteroni</i>	Coma. testosteroni	COMTES
<i>Corynebacterium amycolatum</i>	Cory. amycolatum	CORAMY

Long Name	Short Name	Code
<i>Corynebacterium amycolatum/minuteissimum</i>	Cory. amycolatum/ minuteissimum	CORAMYMIN
<i>Corynebacterium amycolatum/striatum</i>	Cory. amycolatum/striatum	CORAMYSTR
<i>Corynebacterium bovis</i>	Cory. bovis	CORBOV
<i>Corynebacterium diphtheriae</i>	Cory. diphtheriae	CORDIP
<i>Corynebacterium jeikeium</i>	Cory. jeikeium	CORJEI
<i>Corynebacterium kutscheri</i>	Cory. kutscheri	CORKUT
<i>Corynebacterium matruchotii</i>	Cory. matruchotii	CORMAT
<i>Corynebacterium minuteissimum</i>	Cory. minuteissimum	CORMIN
<i>Corynebacterium propinquum</i>	Cory. propinquum	CORPRO
<i>Corynebacterium pseudodiphtheriticum</i>	Cory. pseudodiphth.	CORPSD
<i>Corynebacterium pseudotuberculosis</i>	Cory. pseudotuberc.	CORPST
<i>Corynebacterium renale</i>	Cory. renale	CORREN
<i>Corynebacterium striatum</i>	Cory. striatum	CORSTR
<i>Corynebacterium ulcerans</i>	Cory. ulcerans	CORULC
<i>Corynebacterium urealyticum</i>	Cory. urealyticum	CORURE
<i>Corynebacterium xerosis</i>	Cory. xerosis	CORXER
<i>Cosenzaea myxofaciens</i>	Cosen. myxofaciens	PROTMYX
<i>Cronobacter sakazakii</i> complex	Cronob. sakazakii complex	ENTBSAK
<i>Cryptococcus albidus</i>	Cryp. albidus	CRYALB
<i>Cryptococcus humicola</i>	Cryp. humicola	CRYHUM
<i>Cryptococcus laurentii</i>	Cryp. laurentii	CRYLAU
<i>Cryptococcus luteolus</i>	Cryp. luteolus	CRYLUT
<i>Cryptococcus neoformans</i>	Cryp. neoformans	CRYNEO
<i>Cryptococcus terreus</i>	Cryp. terreus	CRYTER
<i>Cryptococcus uniguttulatus</i>	Cryp. uniguttulatus	CRYUNI
<i>Cupriavidus gilardii</i>	Cup. gilardii	RALGIL

Long Name	Short Name	Code
<i>Cupriavidus pauculus</i>	Cup. pauculus	CDCIVC2
<i>Delftia acidovorans</i>	Delf. acidovorans	COMACI
<i>Dermabacter hominis</i>	Dermab. hominis	DERBHOM
<i>Dermacoccus nishinomiyaensis</i>	Derm. nishinomiyaen.	MICNIS
<i>Edwardsiella hoshinae</i>	Ed. hoshinae	EDWHOS
<i>Edwardsiella ictaluri</i>	Ed. ictaluri	EDWICT
<i>Edwardsiella species</i>	Ed. species	EDWSPE
<i>Edwardsiella tarda</i>	Ed. tarda	EDWTAR
<i>Edwardsiella tarda</i> biogroup 1	Ed. tarda biogr. 1	EDWTAR1
<i>Eikenella corrodens</i>	Eik. corrodens	EIKCOR
<i>Elizabethkingia meningoseptica</i>	Eliz. meningosept.	CHRBMEN
<i>Empedobacter brevis</i>	Emp. brevis	EMPBRE
<i>Enterobacter aerogenes</i>	Enterob. aerogenes	ENTBAER
<i>Enterobacter asburiae</i>	Enterob. asburiae	ENTBASB
<i>Enterobacter cancerogenus</i>	Enterob. cancerogenus	ENTBCAN
<i>Enterobacter cloacae</i>	Enterob. cloacae	ENTBCLO
<i>Enterobacter cloacae</i> ssp <i>dissolvens</i>	Enterob. cloacae ssp dissolven	ENTBDIS
<i>Enterobacter hormaechei</i>	Enterob. hormaechei	ENTBHOR
<i>Enterobacter kobei</i>	Enterob. kobei	ENTBKOB
<i>Enterobacter nimipressuralis</i>	Enterob. nimipressuralis	ENTBNIM
<i>Enterobacter species</i>	Enterob. species	ENTBSPE
<i>Enterococcus asini</i>	Enteroc. asini	ENTCASI
<i>Enterococcus avium</i>	Enteroc. avium	ENTCAVI
<i>Enterococcus casseliflavus</i>	Enteroc. casseliflavus	ENTCCAS
<i>Enterococcus casseliflavus/gallinarum</i>	Enteroc. cassel./gallin.	ENTCCASGAL
<i>Enterococcus cecorum</i>	Enteroc. cecorum	ENTCCEC

Long Name	Short Name	Code
<i>Enterococcus columbae</i>	Enteroc. columbae	ENTCCOL
<i>Enterococcus dispar</i>	Enteroc. dispar	ENTCDIS
<i>Enterococcus durans</i>	Enteroc. durans	ENTCDUR
<i>Enterococcus durans/faecium</i>	Enteroc. durans/faecium	ENTCDURFAI
<i>Enterococcus faecalis</i>	Enteroc. faecalis	ENTCFAA
<i>Enterococcus faecalis/faecium</i>	Enteroc. faecalis/faecium	ENTCFAAFAI
<i>Enterococcus faecium</i>	Enteroc. faecium	ENTCFAI
<i>Enterococcus flavescens</i>	Enteroc. flavescens	ENTCFLA
<i>Enterococcus gallinarum</i>	Enteroc. gallinarum	ENTCGAL
<i>Enterococcus gilvus</i>	Enteroc. gilvus	ENTCGIL
<i>Enterococcus haemoperoxidus</i>	Enteroc. haemoperoxidus	ENTCHAE
<i>Enterococcus hirae</i>	Enteroc. hirae	ENTCHIR
<i>Enterococcus hirae/faecium</i>	Enteroc. hirae/faecium	ENTCHIRFAI
<i>Enterococcus malodoratus</i>	Enteroc. malodoratus	ENTCMAL
<i>Enterococcus moraviensis</i>	Enteroc. moraviensis	ENTCMOR
<i>Enterococcus mundtii</i>	Enteroc. mundtii	ENTCMUN
<i>Enterococcus pallens</i>	Enteroc. pallens	ENTCPAL
<i>Enterococcus pseudoavium</i>	Enteroc. pseudoavium	ENTCPSE
<i>Enterococcus raffinosus</i>	Enteroc. raffinosus	ENTCRAF
<i>Enterococcus raffinosus/avium</i>	Enteroc. raffinosus/avium	ENTCRAFAVI
<i>Enterococcus ratti</i>	Enteroc. ratti	ENTCRAT
<i>Enterococcus saccharolyticus</i>	Enteroc. saccharolyticus	ENTCSAC
<i>Enterococcus species</i>	Enteroc. species	ENTCSPE
<i>Enterococcus sulfureus</i>	Enteroc. sulfureus	ENTCSUL
<i>Erysipelothrix rhusiopathiae</i>	Ery. rhusiopathiae	ERYRHU
<i>Escherichia coli</i>	Esch. coli	ESCCOL

Long Name	Short Name	Code
<i>Escherichia coli</i> serotype O111	Esch. coli O111	ESCCOL0111
<i>Escherichia coli</i> serotype O157	Esch. coli O157	ESCCOL0157
<i>Escherichia fergusonii</i>	Esch. fergusonii	ESCFER
<i>Escherichia hermannii</i>	Esch. hermannii	ESCHER
<i>Escherichia species</i>	Esch. species	ESCSPE
<i>Escherichia vulneris</i>	Esch. vulneris	ESCVUL
<i>Ewingella americana</i>	Ew. americana	EWIAME
<i>Exophiala dermatitidis</i>	Exo. dermatitidis	WANDER
<i>Exophiala species</i>	Exo. species	EXOSPE
<i>Gardnerella vaginalis</i>	Gard. vaginalis	GARVAG
<i>Gemella haemolysans</i>	Gem. haemolysans	GEMHAE
<i>Gemella morbillorum</i>	Gem. morbillorum	GEMMOR
<i>Gemella species</i>	Gem. species	GEMSPE
<i>Geotrichum species</i>	Geo. species	GEOSPE
<i>Globicatella sanguinis</i>	Glob. sanguinis	GLOSAN
Gram-negative rod unidentified enteric	GNR unident. enteric	GNRUNIE
Gram-negative rod unidentified nonfermenter	GNR unident. NF	GNRUNINF
<i>Grimontia hollisae</i>	Grim. hollisae	VIBHOL
<i>Hafnia alvei</i>	Haf. alvei	HAFALV
<i>Hafnia alvei</i> group 1	Haf. alvei gr. 1	HAFALV1
<i>Helcococcus kunzii</i>	Helco. kunzii	HELCKUN
<i>Hortaea werneckii</i>	Hor. werneckii	HORWER
<i>Hyphopichia burtonii</i>	Hyphop. burtonii	PICBUR
<i>Kingella denitrificans</i>	King. denitrificans	KINDEN
<i>Kingella kingae</i>	King. kingae	KINKIN
<i>Klebsiella granulomatis</i>	Kleb. granulomatis	CALYGRA

Long Name	Short Name	Code
<i>Klebsiella oxytoca</i>	Kleb. oxytoca	KLEOXY
<i>Klebsiella ozaenae</i>	Kleb. ozaenae	KLEPNEO
<i>Klebsiella pneumoniae</i>	Kleb. pneumoniae	KLEPNEP
<i>Klebsiella rhinoscleromatis</i>	Kleb. rhinoscleromatis	KLEPNER
<i>Klebsiella species</i>	Kleb. species	KLESPE
<i>Kloeckera species</i>	Kloeck. species	KLOSPE
<i>Kluyvera ascorbata</i>	Kluyvera ascorbata	KLUASC
<i>Kluyvera cryocrescens</i>	Kluyvera cryocres.	KLUCRY
<i>Kluyvera georgiana</i>	Kluyvera georgiana	KLUGEO
<i>Kluyvera intermedia</i>	Kluyvera intermedia	ENTBINT
<i>Kluyvera species</i>	Kluyvera species	KLUSPE
<i>Kocuria kristinae</i>	Koc. kristinae	MICKRI
<i>Kocuria rosea</i>	Koc. rosea	MICROS
<i>Kocuria varians</i>	Koc. varians	MICVAR
<i>Kosakonia cowanii</i>	Kosak. cowanii	ENTBCOW
<i>Kytococcus sedentarius</i>	Kyto. sedentarius	MICSED
<i>Lactococcus garvieae</i>	Lactoc. garvieae	LACCGAR
<i>Lactococcus lactis ssp cremoris</i>	Lactoc. lactis ssp crem.	LACCLACC
<i>Lactococcus lactis ssp hordniae</i>	Lactoc. lactis ssp hord.	LACCLACH
<i>Lactococcus lactis ssp lactis</i>	Lactoc. lactis ssp lactis	LACCLACL
<i>Lactococcus plantarum</i>	Lactoc. plantarum	LACCPLA
<i>Lactococcus raffinolactis</i>	Lactoc. raffinolactis	LACCRAF
<i>Lactococcus species</i>	Lactoc. species	LACCSPE
<i>Leclercia adecarboxylata</i>	Lec. adecarboxylata	LECADE
<i>Leifsonia aquatica</i>	Leif. aquatica	CORAQU
<i>Lelliottia amnigena</i>	Lell. amnigena	ENTBAMN

Long Name	Short Name	Code
<i>Lelliottia amnigena</i> biogroup 1	Lell. amnigena biogr. 1	ENTBAMN1
<i>Lelliottia amnigena</i> biogroup 2	Lell. amnigena biogr. 2	ENTBAMN2
<i>Leminorella grimontii</i>	Lem. grimontii	LEMGR1
<i>Leminorella richardii</i>	Lem. richardii	LEMRIC
<i>Leuconostoc citreum</i>	Leu. citreum	LEUCIT
<i>Leuconostoc lactis</i>	Leu. lactis	LEULAC
<i>Leuconostoc mesenteroides</i> ssp <i>cremoris</i>	Leu. mesenter. ssp crem.	LEUMESC
<i>Leuconostoc mesenteroides</i> ssp <i>mesenteroides</i>	Leu. mesenter. ssp mes.	LEUMESM
<i>Leuconostoc pseudomesenteroides</i>	Leu. pseudomesenter.	LEUPSE
<i>Leuconostoc</i> species	Leu. species	LEUSPE
<i>Listeria grayi</i>	Lis. grayi	LISGRA
<i>Listeria innocua</i>	Lis. innocua	LISINN
<i>Listeria ivanovii</i>	Lis. ivanovii	LISIVA
<i>Listeria monocytogenes</i>	Lis. monocytogenes	LISMON
<i>Listeria monocytogenes/innocua</i>	Lis. monocytogenes/innocua	LISMONINN
<i>Listeria</i> species	Lis. species	LISSPE
<i>Listeria welshimeri</i>	Lis. welshimeri	LISWEL
<i>Lysinibacillus sphaericus</i>	Lysini. sphaericus	BACISPH
<i>Macrococcus caseolyticus</i>	Macro. caseolyticus	STACAS
<i>Magnusiomyces capitatus</i>	Magnus. capitatus	BLACAP
<i>Malassezia furfur</i> complex	Mal. furfur complex	MALFURX
<i>Malassezia pachydermatis</i>	Mal. pachydermatis	MALPAC
<i>Malassezia sympodialis</i>	Mal. sympodialis	MALSYM
<i>Mannheimia haemolytica</i>	Mann. haemolytica	PASHAE
<i>Methylobacterium extorquens</i>	Methylob. extorquens	METEXT
<i>Micrococcus luteus</i>	Microc. luteus	MICLUT

Long Name	Short Name	Code
<i>Micrococcus lylae</i>	Microc. lylae	MICLYL
<i>Millerozyma farinosa</i>	Mill. farinosa	PICFAR
<i>Moellerella wisconsensis</i>	Moel. wisconsensis	MOEWIS
<i>Moraxella (Branhamella) catarrhalis</i>	Morax. (Bran.) cat.	MORABRACAT
<i>Moraxella species</i>	Morax. species	MORASPE
<i>Morganella morganii</i>	Morg. morganii	MORGMOR
<i>Morganella morganii ssp morganii</i>	Morg. morg. ssp morg.	MORGMORM
<i>Morganella morganii ssp morganii biogroup 1</i>	Morg. morg. ssp morg. biog. 1	MORGMORM1
<i>Morganella morganii ssp sibonii</i>	Morg. morg. ssp sibonii	MORGMORS
<i>Morganella species</i>	Morg. species	MORGspe
<i>Myroides odoratimimus</i>	Myr. odoratimimus	MYRODI
<i>Myroides odoratus</i>	Myr. odoratus	MYRODA
<i>Myroides odoratus/odoratimimus</i>	Myr. odoratus/odoratimimus	MYRODAODI
<i>Myroides species</i>	Myr. species	MYRSPE
<i>Neisseria animaloris</i>	Nei. animaloris	CDCEF4a
<i>Neisseria zoodegmatidis</i>	Nei. zoodegmatidis	CDCEF4b
<i>Ochrobactrum anthropi</i>	Och. anthropi	OCHANT
<i>Oligella ureolytica</i>	Olig. ureolytica	OLIURO
<i>Oligella urethralis</i>	Olig. urethralis	OLIURT
<i>Paenibacillus alvei</i>	Paen. alvei	PAEALV
<i>Paenibacillus macerans</i>	Paen. macerans	PAEMAC
<i>Pantoea agglomerans</i>	Pan. agglomerans	PANAGG
<i>Pantoea ananatis</i>	Pan. ananatis	ERWANA
<i>Pantoea dispersa</i>	Pan. dispersa	PANDIS
<i>Pantoea species</i>	Pan. species	PANSPE
<i>Pantoea stewartii</i>	Pan. stewartii	PANSTE

Long Name	Short Name	Code
<i>Pantoea stewartii</i> ssp <i>indologenes</i>	Pan. ste. ssp indologenes	PANSTEI
<i>Pantoea stewartii</i> ssp <i>stewartii</i>	Pan. ste. ssp stewartii	PANSTES
<i>Paracoccus yeei</i>	Parac. yeei	CDCEO2
<i>Pasteurella aerogenes</i>	Past. aerogenes	PASAER
<i>Pasteurella multocida</i>	Past. multocida	PASMUL
<i>Pasteurella pneumotropica</i>	Past. pneumotropica	PASPNE
<i>Pediococcus acidilactici</i>	Ped. acidilactici	PEDACI
<i>Pediococcus damnosus</i>	Ped. damnosus	PEDDAM
<i>Pediococcus dextrinicus</i>	Ped. dextrinicus	PEDDEX
<i>Pediococcus parvulus</i>	Ped. parvulus	PEDPAR
<i>Pediococcus pentosaceus</i>	Ped. pentosaceus	PEDPEN
<i>Pediococcus species</i>	Ped. species	PEDSPE
<i>Photobacterium damsela</i>	Photob. damsela	PHOBDAM
<i>Plesiomonas shigelloides</i>	Ples. shigelloides	PLESHI
<i>Pluralibacter gergoviae</i>	Plural. gergoviae	ENTBGER
<i>Pragia fontium</i>	Prag. fontium	PRAFON
<i>Proteus hauseri</i>	Prot. hauseri	PROTHAU
<i>Proteus mirabilis</i>	Prot. mirabilis	PROTMIR
<i>Proteus penneri</i>	Prot. penneri	PROTPEN
<i>Proteus species</i>	Prot. species	PROTSPE
<i>Proteus vulgaris</i>	Prot. vulgaris	PROTVUL
<i>Proteus vulgaris/penneri</i>	Prot. vulgaris/penneri	PROTVULPEN
<i>Prototheca wickerhamii</i>	Protot. wick.	PROHWIC
<i>Prototheca zopfii</i>	Protot. zopfii	PROHZOP
<i>Providencia alcalifaciens</i>	Prov. alcalifaciens	PROVALC
<i>Providencia heimbachae</i>	Prov. heimbachae	PROVHEI

Long Name	Short Name	Code
<i>Providencia rettgeri</i>	Prov. rettgeri	PROVRET
<i>Providencia rustigianii</i>	Prov. rustigianii	PROVRUS
<i>Providencia species</i>	Prov. species	PROVSPE
<i>Providencia stuartii</i>	Prov. stuartii	PROVSTU
<i>Pseudomonas aeruginosa</i>	Pseud. aeruginosa	PSEAER
<i>Pseudomonas alcaligenes</i>	Pseud. alcaligenes	PSEALC
<i>Pseudomonas fluorescens</i>	Pseud. fluorescens	PSEFLU
<i>Pseudomonas fluorescens/putida</i>	Pseud. fluorescens/putida	PSEFLUPUT
<i>Pseudomonas luteola</i>	Pseud. luteola	CHRMLUT
<i>Pseudomonas mendocina</i>	Pseud. mendocina	PSEMEN
<i>Pseudomonas monteilii</i>	Pseud. monteilii	PSEMON
<i>Pseudomonas oryzihabitans</i>	Pseud. oryzihabitans	FLAIORY
<i>Pseudomonas pertucinogena</i>	Pseud. pertucinogena	PSEPER
<i>Pseudomonas pseudoalcaligenes</i>	Pseud. pseudoalcaligenes	PSEPSE
<i>Pseudomonas pseudoalcaligenes ssp pseudoalcaligenes</i>	Pseud. pseudoal. ssp pseud.	PSEPSEP
<i>Pseudomonas putida</i>	Pseud. putida	PSEPUT
<i>Pseudomonas species</i>	Pseud. species	PSESPE
<i>Pseudomonas stutzeri</i>	Pseud. stutzeri	PSESTU
<i>Pseudomonas veronii</i>	Pseud. veronii	PSEVER
<i>Rahnella aquatilis</i>	Rah. aquatilis	RAHAQU
<i>Ralstonia pickettii</i>	Ral. pickettii	BURPIC
<i>Ralstonia solanacearum</i>	Ral. solanacearum	BURSOL
<i>Ralstonia species</i>	Ral. species	RALSPE
<i>Raoultella ornithinolytica</i>	Rao. ornithinolytica	KLEORN
<i>Raoultella planticola</i>	Rao. planticola	KLEPLA
<i>Raoultella species</i>	Rao. species	RAOSPE

Long Name	Short Name	Code
<i>Raoultella terrigena</i>	Rao. terrigena	KLETER
<i>Rhizobium radiobacter</i>	Rhizob. radiobacter	AGRRAD
<i>Rhodococcus equi</i>	Rhodoc. equi	RHOCEQU
<i>Rhodotorula glutinis</i>	Rhodot. glutinis	RHOTGLU
<i>Rhodotorula minuta</i>	Rhodot. minuta	RHOTMIN
<i>Rhodotorula mucilaginosa</i> var <i>mucilaginosa</i>	Rhodot. mucilag. var mucilag.	RHOTMUCM
<i>Rothia dentocariosa</i>	Roth. dentocariosa	ROTDEN
<i>Rothia mucilaginosa</i>	Roth. mucilaginosa	STOMUC
<i>Saccharomyces cerevisiae</i>	Sac. cerevisiae	SACCER
<i>Salmonella enterica</i> ssp <i>arizonae</i>	Salm. enterica ssp arizonae	SALCHOA
<i>Salmonella enterica</i> ssp <i>diarizonae</i>	Salm. enterica ssp diarizonae	SALCHOD
<i>Salmonella enterica</i> ssp <i>enterica</i> serovar <i>Choleraesuis</i>	Salm. enterica sv Choleraesuis	SALCHOC
<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Gallinarum</i> bv <i>Gallinarum</i>	Salm. Gallinarum	SALGAL
<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Gallinarum</i> bv <i>Pullorum</i>	Salm. Pullorum	SALPUL
<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Paratyphi</i> A	Salm. Paratyphi A	SALPARA
<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Typhi</i>	Salm. Typhi	SALTYP
<i>Salmonella enterica</i> ssp <i>houtenae</i>	Salm. enterica ssp houtenae	SALCHOH
<i>Salmonella enterica</i> ssp <i>indica</i>	Salm. enterica ssp indica	SALCHOI
<i>Salmonella enterica</i> ssp <i>salamae</i>	Salm. enterica ssp salamae	SALCHOS
<i>Salmonella</i> species	Salm. species	SALSPE
<i>Serratia entomophila</i>	Ser. entomophila	SERENT
<i>Serratia ficaria</i>	Ser. ficaria	SERFIC
<i>Serratia fonticola</i>	Ser. fonticola	SERFON
<i>Serratia grimesii</i>	Ser. grimesii	SERGRI
<i>Serratia liquefaciens</i>	Ser. liquefaciens	SERLIQ

Long Name	Short Name	Code
<i>Serratia marcescens</i>	Ser. marcescens	SERMAR
<i>Serratia odorifera</i>	Ser. odorifera	SERODO
<i>Serratia odorifera</i> biogroup 1	Ser. odorifera 1	SERODO1
<i>Serratia odorifera</i> biogroup 2	Ser. odorifera 2	SERODO2
<i>Serratia plymuthica</i>	Ser. plymuthica	SERPLY
<i>Serratia proteamaculans</i> ssp <i>proteamaculans</i>	Ser. proteamac. ssp proteam.	SERPROP
<i>Serratia proteamaculans</i> ssp <i>quinovora</i>	Ser. proteamac. ssp quino.	SERPROQ
<i>Serratia rubidaea</i>	Ser. rubidaea	SERRUB
<i>Serratia</i> species	Ser. species	SERSPE
<i>Shewanella algae</i>	Shew. algae	SHEALG
<i>Shewanella putrefaciens</i>	Shew. putrefaciens	SHEPUT
<i>Shewanella</i> species	Shew. species	SHESPE
<i>Shigella boydii</i>	Shig. boydii	SHIBOY
<i>Shigella dysenteriae</i>	Shig. dysenteriae	SHIDYS
<i>Shigella flexneri</i>	Shig. flexneri	SHIFLE
<i>Shigella sonnei</i>	Shig. sonnei	SHISON
<i>Shigella</i> species	Shig. species	SHISPE
<i>Shimwellia blattae</i>	Shim. blattae	ESCBLA
<i>Sphingobacterium multivorum</i>	Sphingob. multivorum	SPHBMUL
<i>Sphingobacterium multivorum/thalpophilum</i>	Sphb. multivorum/thalpophilum	SPHBMULTHA
<i>Sphingobacterium</i> species	Sphingob. species	SPHBSPE
<i>Sphingobacterium spiritivorum</i>	Sphingob. spiritivorum	SPHBSPi
<i>Sphingobacterium thalpophilum</i>	Sphingob. thalpophilum	SPHBTHA
<i>Sphingomonas paucimobilis</i>	Sphingom. paucimobilis	SPHMPAU
<i>Sporobolomyces salmonicolor</i>	Sporobol. salmonicolor	SPOBSAL
<i>Staphylococcus arlettae</i>	Staph. arlettae	STAARL

Long Name	Short Name	Code
<i>Staphylococcus aureus</i>	Staph. aureus	STAAUE
<i>Staphylococcus aureus ssp anaerobius</i>	Staph. aureus ssp anaerob.	STAAUEAN
<i>Staphylococcus aureus ssp aureus</i>	Staph. aureus ssp aureus	STAAUEAU
<i>Staphylococcus auricularis</i>	Staph. auricularis	STAAUI
<i>Staphylococcus capitis</i>	Staph. capitis	STACAI
<i>Staphylococcus capitis ssp capitis</i>	Staph. capitis ssp capitis	STACAIC
<i>Staphylococcus capitis ssp urealyticus</i>	Staph. capitis ssp urealyt.	STACAIU
<i>Staphylococcus caprae</i>	Staph. caprae	STACAP
<i>Staphylococcus carnosus</i>	Staph. carnosus	STACAR
<i>Staphylococcus carnosus ssp carnosus</i>	Staph. carn. ssp carn.	STACARC
<i>Staphylococcus carnosus ssp utilis</i>	Staph. carn. ssp utilis	STACARU
<i>Staphylococcus chromogenes</i>	Staph. chromogenes	STACHR
<i>Staphylococcus chromogenes/hyicus</i>	Staph. chromogenes/hyicus	STACHRHYI
<i>Staphylococcus coagulase-negative</i>	Staph. coag. neg.	STACNEG
<i>Staphylococcus coagulase-positive</i>	Staph. coag. pos.	STACPOS
<i>Staphylococcus cohnii</i>	Staph. cohnii	STACOH
<i>Staphylococcus cohnii ssp cohnii</i>	Staph. cohnii ssp cohnii	STACOHC
<i>Staphylococcus cohnii ssp urealyticum</i>	Staph. cohnii ssp urealyt.	STACOHU
<i>Staphylococcus condimentii</i>	Staph. condimentii	STACON
<i>Staphylococcus delphini</i>	Staph. delphini	STADEL
<i>Staphylococcus epidermidis</i>	Staph. epidermidis	STAEPI
<i>Staphylococcus equorum</i>	Staph. equorum	STAEQU
<i>Staphylococcus felis</i>	Staph. felis	STAFEL
<i>Staphylococcus fleurettii</i>	Staph. fleurettii	STAFLE
<i>Staphylococcus gallinarum</i>	Staph. gallinarum	STAGAL
<i>Staphylococcus haemolyticus</i>	Staph. haemolyticus	STAHAE

Long Name	Short Name	Code
<i>Staphylococcus haemolyticus/lugdunensis</i>	Staph. haemol./lugdun.	STAHAE LUG
<i>Staphylococcus hominis</i>	Staph. hominis	STAHOM
<i>Staphylococcus hominis ssp hominis</i>	Staph. hom. ssp hom.	STAHOMH
<i>Staphylococcus hominis ssp novobiosepticus</i>	Staph. hom. ssp novo.	STAHOMN
<i>Staphylococcus hyicus</i>	Staph. hyicus	STAHYI
<i>Staphylococcus intermedius</i>	Staph. intermedius	STAI NT
<i>Staphylococcus kloosii</i>	Staph. kloosii	STAKLO
<i>Staphylococcus lentus</i>	Staph. lentus	STALEN
<i>Staphylococcus lugdunensis</i>	Staph. lugdunensis	STALUG
<i>Staphylococcus lutrae</i>	Staph. lutrae	STALUT
<i>Staphylococcus muscae</i>	Staph. muscae	STAMUS
<i>Staphylococcus pasteurii</i>	Staph. pasteurii	STAPAS
<i>Staphylococcus pettenkoferi</i>	Staph. pettenkoferi	STAPET
<i>Staphylococcus piscifermentans</i>	Staph. piscifermentans	STAPIS
<i>Staphylococcus pulvereri</i>	Staph. pulvereri	STAPUL
<i>Staphylococcus saccharolyticus</i>	Staph. saccharolyticus	STASAC
<i>Staphylococcus saprophyticus</i>	Staph. saprophyticus	STASAP
<i>Staphylococcus saprophyticus ssp bovis</i>	Staph. sap. ssp bovis	STASAPB
<i>Staphylococcus saprophyticus ssp saprophyticus</i>	Staph. sap. ssp saprophyticus	STASAPS
<i>Staphylococcus schleiferi</i>	Staph. schleiferi	STASCH
<i>Staphylococcus schleiferi ssp coagulans</i>	Staph. schleiferi ssp coagul.	STASCHC
<i>Staphylococcus schleiferi ssp schleiferi</i>	Staph. schleiferi ssp schleif.	STASCHS
<i>Staphylococcus sciuri</i>	Staph. sciuri	STASCI
<i>Staphylococcus sciuri ssp carnaticus</i>	Staph. sciuri ssp carnaticus	STASCIC
<i>Staphylococcus sciuri ssp rodentium</i>	Staph. sciuri ssp rodentium	STASCIR
<i>Staphylococcus sciuri ssp sciuri</i>	Staph. sciuri ssp sciuri	STASCIS

Long Name	Short Name	Code
<i>Staphylococcus simulans</i>	Staph. simulans	STASIM
<i>Staphylococcus species</i>	Staph. species	STASPE
<i>Staphylococcus succinus</i>	Staph. succinus	STASUC
<i>Staphylococcus succinus ssp casei</i>	Staph. suc. ssp casei	STASUCCA
<i>Staphylococcus succinus ssp succinus</i>	Staph. suc. ssp succinus	STASUCSU
<i>Staphylococcus vitulinus</i>	Staph. vitulinus	STAVIT
<i>Staphylococcus warneri</i>	Staph. warneri	STAWAR
<i>Staphylococcus warneri/pasteuri</i>	Staph. warneri/pasteuri	STAWARPAS
<i>Staphylococcus xylosus</i>	Staph. xylosus	STAXYL
<i>Stenotrophomonas africana</i>	Sten. africana	STEAFR
<i>Stenotrophomonas maltophilia</i>	Sten. maltophilia	STEMAL
<i>Streptococcus acidominimus</i>	Strep. acidominimus	STRACI
<i>Streptococcus agalactiae</i> (Strep. group B)	Strep. agalactiae (Str. gr. B)	STRAGA
<i>Streptococcus alactolyticus</i>	Strep. alactolyticus	STRALA
<i>Streptococcus alpha-hemolytic</i>	Strep. alpha-hemolytic	STRAHE
<i>Streptococcus anginosus</i>	Strep. anginosus	STRANG
<i>Streptococcus anginosus</i> (previously milleri) group	Strep. anginosus (milleri) gr.	STRANGGR
<i>Streptococcus beta-hemolytic</i> ACG (large colony)	Strep. beta-hemo ACG (lg col)	STRBHE
<i>Streptococcus bovis</i> (Strep. group D)	Strep. bovis (Str. gr. D)	STRBOV
<i>Streptococcus bovis</i> I (Strep. group D)	Strep. bovis I (Str. gr. D)	STRBOVI
<i>Streptococcus bovis</i> II (Strep. group D)	Strep. bovis II (Str. gr. D)	STRBOVII
<i>Streptococcus canis</i>	Strep. canis	STRCAN
<i>Streptococcus constellatus</i>	Strep. constellatus	STRCON
<i>Streptococcus constellatus ssp constellatus</i>	Strep. con ssp constellatus	STRCONCO
<i>Streptococcus constellatus ssp pharyngis</i>	Strep. con ssp pharyngis	STRCONPH
<i>Streptococcus criceti</i>	Strep. criceti	STRCRC

Long Name	Short Name	Code
<i>Streptococcus cristatus</i>	Strep. cristatus	STRCRS
<i>Streptococcus downei</i>	Strep. downei	STRDOW
<i>Streptococcus dysgalactiae</i>	Strep. dysgalactiae	STRDYS
<i>Streptococcus dysgalactiae</i> ssp <i>dysgalactiae</i>	Strep. dysgal. ssp dysgal.	STRDYSDY
<i>Streptococcus dysgalactiae</i> ssp <i>equisimilis</i>	Strep. dysgal. ssp equis.	STRDYSEM
<i>Streptococcus dysgalactiae/canis</i>	Strep. dysgal./canis	STRDYSCAN
<i>Streptococcus equi</i>	Strep. equi	STREQU
<i>Streptococcus equi</i> ssp <i>equi</i>	Strep. equi ssp equi	STREQUE
<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>	Strep. equi ssp zooepid.	STREQUZ
<i>Streptococcus equinus</i>	Strep. equinus	STREQN
<i>Streptococcus ferus</i>	Strep. ferus	STRFER
<i>Streptococcus gallolyticus</i> ssp <i>macedonicus</i>	Strep. gallolyti. ssp macedon.	STRMAC
<i>Streptococcus gordonii</i>	Strep. gordonii	STRGOR
<i>Streptococcus</i> group A (small colony)	Strep. group A (sm col)	STRGRAS
<i>Streptococcus</i> group A (<i>Strep. pyogenes</i>)	Strep. group A (Str. pyogenes)	STRGRA
<i>Streptococcus</i> group B (<i>Strep. agalactiae</i>)	Strep. group B (Str. agalact.)	STRGRB
<i>Streptococcus</i> group C (large colony)	Strep. group C (lg col)	STRGRC
<i>Streptococcus</i> group C (small colony)	Strep. group C (sm col)	STRGRCS
<i>Streptococcus</i> group C/G (large colony)	Strep. group C/G (lg col)	STRGRCG
<i>Streptococcus</i> group C/G (small colony)	Strep. group C/G (sm col)	STRGRCGS
<i>Streptococcus</i> group CFG (small colony)	Strep. group CFG (sm col)	STRGRCFG
<i>Streptococcus</i> group D (non-enterococcus)	Strep. group D (non-enteroc.)	STRGRDNE
<i>Streptococcus</i> group E	Strep. group E	STRGRE
<i>Streptococcus</i> group F	Strep. group F	STRGRF
<i>Streptococcus</i> group G (large colony)	Strep. group G (lg col)	STRGRG
<i>Streptococcus</i> group G (small colony)	Strep. group G (sm col)	STRGRGS

Long Name	Short Name	Code
<i>Streptococcus</i> group L	Strep. group L	STRGRL
<i>Streptococcus hyointestinalis</i>	Strep. hyointestinalis	STRHYO
<i>Streptococcus infantarius</i>	Strep. infantarius	STRINA
<i>Streptococcus infantarius ssp coli</i>	Strep. infa ssp coli	STRINACO
<i>Streptococcus infantarius ssp infantarius</i>	Strep. infa ssp infantarius	STRINAIN
<i>Streptococcus infantis</i>	Strep. infantis	STRINF
<i>Streptococcus iniae</i>	Strep. iniae	STRINI
<i>Streptococcus intermedius</i>	Strep. intermedius	STRINR
<i>Streptococcus milleri</i> group	Strep. milleri gr.	STRMILGR
<i>Streptococcus mitis</i>	Strep. mitis	STRMIT
<i>Streptococcus mitis</i> group	Strep. mitis gr.	STRMITGR
<i>Streptococcus mitis/oralis</i>	Strep. mitis/oralis	STRMITORA
<i>Streptococcus mitis/pneumoniae</i>	Strep. mitis/pneum.	STRMITPNE
<i>Streptococcus mutans</i>	Strep. mutans	STRMUT
<i>Streptococcus mutans</i> group	Strep. mutans gr.	STRMUTGR
<i>Streptococcus oralis</i>	Strep. oralis	STRORA
<i>Streptococcus parasanguinis</i>	Strep. parasanguinis	STRPAR
<i>Streptococcus peroris</i>	Strep. peroris	STRPER
<i>Streptococcus pneumoniae</i>	Strep. pneumoniae	STRPNE
<i>Streptococcus porcinus</i>	Strep. porcinus	STRPOR
<i>Streptococcus pyogenes</i> (Strep. group A)	Strep. pyogenes (Str. gr. A)	STRPYO
<i>Streptococcus ratti</i>	Strep. ratti	STRRAT
<i>Streptococcus salivarius</i>	Strep. salivarius	STRSAL
<i>Streptococcus salivarius</i> group	Strep. salivarius gr.	STRSALGR
<i>Streptococcus salivarius ssp thermophilus</i>	Strep. salivar. ssp thermoph.	STRTHE
<i>Streptococcus sanguinis</i>	Strep. sanguinis	STRSAN

Long Name	Short Name	Code
<i>Streptococcus sanguinis</i> group	Strep. sanguinis gr.	STRSANGR
<i>Streptococcus sobrinus</i>	Strep. sobrinus	STRSOB
<i>Streptococcus species</i>	Strep. species	STRSPE
<i>Streptococcus suis</i>	Strep. suis	STRSUI
<i>Streptococcus uberis</i>	Strep. uberis	STRUBE
<i>Streptococcus vestibularis</i>	Strep. vestibularis	STRVES
<i>Streptococcus viridans</i> beta-hemolytic (small colony)	Strep. vir. beta-hemo (sm col)	STRBHES
<i>Streptococcus viridans</i> group	Strep. viridans gr.	STRVIRGR
<i>Suttonella indologenes</i>	Sutto. indologenes	SUTIND
<i>Tatumella ptyseos</i>	Tat. ptyseos	TATPTY
<i>Trichosporon asahii</i>	Tric. asahii	TRIASA
<i>Trichosporon inkin</i>	Tric. inkin	TRIINK
<i>Trichosporon loubieri</i>	Tric. loubieri	TRILOU
<i>Trichosporon mucoides</i>	Tric. mucoides	TRIMUC
<i>Trichosporon ovoides</i>	Tric. ovoides	TRIOVO
<i>Trueperella pyogenes</i>	True. pyogenes	ACTMPYO
<i>Vibrio alginolyticus</i>	Vib. alginolyticus	VIBALG
<i>Vibrio cholerae</i>	Vib. cholerae	VIBCHO
<i>Vibrio fluvialis</i>	Vib. fluvialis	VIBFLU
<i>Vibrio metschnikovii</i>	Vib. metschnikovii	VIBMET
<i>Vibrio mimicus</i>	Vib. mimicus	VIBMIM
<i>Vibrio parahaemolyticus</i>	Vib. parahaemolyticus	VIBPAR
<i>Vibrio vulnificus</i>	Vib. vulnificus	VIBVUL
<i>Weeksella virosa</i>	Week. virosa	WEEVIR
<i>Yersinia aldovae</i>	Yer. aldovae	YERALD
<i>Yersinia bercovieri</i>	Yer. bercovieri	YERBER

Long Name	Short Name	Code
<i>Yersinia enterocolitica</i>	Yer. enterocolitica	YERENT
<i>Yersinia enterocolitica</i> group	Yer. enterocolitica gr.	YERENTGR
<i>Yersinia frederiksenii</i>	Yer. frederiksenii	YERFRE
<i>Yersinia intermedia</i>	Yer. intermedia	YERINT
<i>Yersinia kristensenii</i>	Yer. kristensenii	YERKRI
<i>Yersinia mollaretii</i>	Yer. mollaretii	YERMOL
<i>Yersinia pseudotuberculosis</i>	Yer. pseudotuberculosis	YERPSE
<i>Yersinia rohdei</i>	Yer. rohdei	YERROH
<i>Yersinia ruckeri</i>	Yer. ruckeri	YERRUC
<i>Yersinia species</i>	Yer. species	YERSPE
<i>Yokenella regensburgei</i>	Yok. regensburgei	YOKREG
<i>Zygosaccharomyces bailii</i>	Zyg. bailii	ZYGBAI

Panel Information

7.4 List of Antimicrobial Agents in BD Phoenix Panels

7.4.1 Gram Negative

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Aminoglycoside	Amikacin	AN	0.5-64	0.5-64	N/A
Beta-lactam	Amoxicillin	AMX	0.5-32	0.25-32	0.125-32
Beta-lactam	Amoxicillin-Clavulanate (f)	AXC	0.5/2-32/2	0.25/2-32/2	N/A
Beta-lactam	Amoxicillin-Clavulanate	AMC	0.5/0.25-32/16	0.25/0.125-32/16	0.125/0.06-32/16
Beta-lactam	Ampicillin	AM	0.5-32	0.0625-32	0.0625-32
Beta-lactam	Ampicillin-Sulbactam (f)	SXA	0.5/8-32/8	0.5/8-32/8	N/A
Beta-lactam	Ampicillin-Sulbactam	SAM	1/0.5-32/16	0.5/0.25-32/16	N/A
Aminoglycoside	Arbekacin	ARB	0.25-16	0.25-16	N/A
Macrolides Lincosmides Streptogramins	Azithromycin	AZM	N/A	0.06-8	N/A
Beta-lactam	Aztreonam	ATM	0.5-64	N/A	N/A
Beta-lactam	Cefaclor	CEC	N/A	0.5-32	N/A
Beta-lactam	Cefazolin	CZ	0.5-32	0.5-32	N/A
Beta-lactam	Cefdinir	CDR	0.125-4	0.125-4	N/A
Beta-lactam	Cefditoren	CDN	0.125-8	0.125-8	N/A
Beta-lactam	Cefepime	FEP	0.125-64	1-64	0.0625-4
Beta-lactam	Cefetamet-pivoxil	CAT	0.25-16	N/A	N/A
Beta-lactam	Cefixime	CFM	0.125-8	N/A	N/A

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Beta-lactam	Cefmetazole	CMZ	2-64	1-64	N/A
Beta-lactam	Cefoperazone	CFP	0.5-64	1-64	N/A
Beta-lactam	Cefoperazone-Sulbactam	SCP	0.5/8-64/8	N/A	N/A
Beta-lactam	Cefotaxime	CTX	0.5-64	1-64	0.0625-4
Beta-lactam	Cefotetan	CTT	2-64	1-64	N/A
Beta-lactam	Cefotiam	CFT	0.5-64	0.5-64	N/A
Beta-lactam	Cefoxitin	FOX	0.5-64	1-64	N/A
Beta-lactam	Cefozopran	CFZ	N/A	1-64	N/A
Beta-lactam	Cefpirome	CPO	0.5-64	0.5-64	N/A
Beta-lactam	Cefpodoxime	CPD	0.125-8	0.5-8	N/A
Beta-lactam	Cefsulodin	CFS	1-64	N/A	N/A
Beta-lactam	Ceftaroline	CPT	N/A	0.0625-4	N/A
Beta-lactam	Ceftazidime	CAZ	0.5-64	1-64	N/A
Beta-lactam	Ceftibuten	CTB	0.5-32	N/A	N/A
Beta-lactam	Ceftizoxime	ZOX	0.5-64	1-64	N/A
Beta-lactam	Ceftriaxone	CRO	0.5-64	1-64	0.0625-4
Beta-lactam	Cefuroxime	CXM	1-64	1-64	0.125-4
Beta-lactam	Cephalexin	CN	1-64	0.5-64	N/A
Beta-lactam	Cephalothin	CF	1-64	0.5-64	N/A
Phenicol	Chloramphenicol	C	1-32	1-32	1-32
Quinolone	Ciprofloxacin	CIP	0.125-4	0.125-4	N/A
Macrolides Lincosamides Streptogramins	Clarithromycin	CLR	N/A	0.06-8	N/A
Macrolide Lincosamide	Clindamycin	CC	N/A	0.12-8	0.03-4
Cyclic peptide	Colistin	CL	0.5-4	N/A	N/A
Cyclic lipopeptide	Daptomycin	DAP	N/A	0.125-32	0.03-16
Tetracycline	Doxycycline	D	N/A	0.25-16	N/A

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Beta-lactam	Ertapenem	ETP	0.0625-4	0.25-32	0.0625-4
Macrolide Lincosamide	Erythromycin	E	N/A	0.125-8	0.015-4
Fosfomycin	Fosfomycin	FF	16-256	8-256	N/A
Fusidane	Fusidic Acid	FA	N/A	0.5-32	N/A
Quinolone	Garenoxacin	GRN	0.125-16	0.25-8	0.03125-4
Quinolone	Gatifloxacin	GAT	0.25-8	0.25-8	0.0625-8
Quinolone	Gemifloxacin	GEM	0.125-8	0.125-2	0.0625-8
Aminoglycoside	Gentamicin	GM	0.5-16	0.5-16	N/A
Aminoglycoside	Gentamicin-Syn	GMS	N/A	500	250-1000
Aminoglycoside	Gentamicin-Syn (SFM)	GMF	N/A	250	N/A
Beta-lactam	Imipenem	IPM	0.0625-16	0.5-16	0.015-4
Aminoglycoside	Isepamycin	ISP	0.5-32	N/A	N/A
Aminoglycoside	Kanamycin	K	0.5-64	0.5-64	N/A
Aminoglycoside	Kanamycin-Syn	KS	N/A	250	250-1000
Quinolone	Levofloxacin	LVX	0.25-8	0.25-8	0.25-16
Macrolides Lincosamides Streptogramin	Lincomycin	L	N/A	0.5-16	N/A
Oxazolidinone	Linezolid	LZD	N/A	0.25-32	0.25-16
Quinolone	Lomefloxacin	LOM	0.25-8	0.25-8	N/A
Beta-lactam	Mecillinam	MEC	0.5-32	N/A	N/A
Beta-lactam	Meropenem	MEM	0.125-32	0.5-16	0.03125-2
Tetracycline	Minocycline	MI	0.5-16	1-32	N/A
Beta-lactam	Moxalactam	MOX	1-64	1-64	N/A
Quinolone	Moxifloxacin	MXF	0.125-8	0.25-8	0.0625-8
Pseudomonic acid	Mupirocin	MUP	N/A	0.0625-8	N/A
Pseudomonic acid	Mupirocin High level	MUH	N/A	256	N/A

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep
Quinolone	Nalidixic Acid	NA	1-32	N/A	N/A
Aminoglycoside	Netilmicin	NET	0.5-32	0.5-32	N/A
Nitrofuran	Nitrofurantoin	FM	8-512	16-512	N/A
Quinolone	Norfloxacin	NOR	0.25-16	0.25-16	N/A
Quinolone	Ofloxacin	OFX	0.25-8	0.25-8	0.5-16
Beta-lactam	Oxacillin	OX	N/A	0.06-4	N/A
Quinolone	Pefloxacin	PEF	0.25-8	0.25-8	N/A
Beta-lactam	Penicillin G	P	N/A	0.0625-32	0.03125-8
Beta-lactam	Piperacillin	PIP	0.5-128	1-128	N/A
Beta-lactam	Piperacillin-Tazobactam	TZP	0.5/4-128/4	1/4-128/4	N/A
Macrolides Lincosamides Streptogramin	Pristinamycin	PR	N/A	0.25-4	0.0625-4
Macrolides Lincosamides Streptogramin	Quinupristin-Dalfopristin	SYN	N/A	0.5-4	0.125-8
Rifamycin	Rifampin	RA	N/A	0.25-32	N/A
Aminoglycoside	Streptomycin-Syn	STS	N/A	1000	250-1000
Aminoglycoside	Streptomycin-Syn (SFM)	STF	N/A	250	N/A
Glycopeptide	Teicoplanin	TEC	N/A	0.5-32	1-32
Ketolide	Telithromycin	TEL	N/A	0.03125-8	0.0625-4
Beta-lactam	Temocillin	TEM	2 – 32	N/A	N/A
Tetracyclinee	Tetracycline	TE	0.5-16	0.5-16	0.0625-16
Beta-lactam	Ticarcillin	TIC	1-128	1-128	N/A
Beta-lactam	Ticarcillin-Clavulanate	TIM	1/2-128/2	1/2-128/2	N/A
Glycylcycline	Tigecycline	TGC	0.25-16	0.03125-4	N/A
Aminoglycoside	Tobramycin	NN	0.12-16	1 -16	N/A

Antimicrobial			Available Concentrations		
			Gram Negative	Gram Positive	Strep

Folate Antagonist	Trimeth-Sulfa (DIN)	STG	0.4/7.6-12.8/243.2	0.4/7.6-12.8/243.2	N/A
Folate Antagonist	Trimethoprim	TMP	0.5-16	0.5-16	N/A
Folate Antagonist	Trimethoprim-Sulfamethoxazole	SXT	0.5/9.5-16/304	0.5/9.5-16/304	0.06/1.19-16/304
Glycopeptide	Vancomycin	VA	N/A	0.5-32	0.0625-32

BOLD = Different range for Gram Negative and Gram Positive

NOTE

MIC dilutions appearing in this manual are actual serial 2-fold dilution concentrations. MIC values appearing on reports may be rounded.

In general, the BD Phoenix M50 instrument provides a MIC for all organisms at any of the concentrations defined on a specific panel. For certain drug and organism combinations in the following tables, a specific minimum/maximum MIC is reported even if there is a lower/higher concentration on the panel.

Antimicrobial	Organism	Min MIC (µg/mL)	Max MIC (µg/mL)
Gram Negative			
Amikacin	<i>Morganella morganii</i>	2	
Amoxicillin-Clavulanic Acid (fixed clav. acid conc.)	<i>Proteus mirabilis</i>	1/2	
Amoxicillin-Clavulanic Acid (variable clav. acid conc.)	<i>Proteus mirabilis</i>	2/1	
Arbekacin	<i>Enterobacter aerogenes</i>	1	
	<i>Morganella morganii</i>	1	
	<i>Providencia</i> genus	1	
Cefdinir	<i>Proteus mirabilis</i>	0.25	

Antimicrobial	Organism	Min MIC (µg/mL)	Max MIC (µg/mL)
Cefixime	<i>Proteus mirabilis</i>	0.5	
	<i>Providencia</i> genus	0.5	
Cefotiam	<i>Citrobacter freundii</i>	2	
	<i>Proteus mirabilis</i>	4	
	<i>Providencia</i> genus	4	
Fosfomycin w/G6P	<i>Pseudomonas aeruginosa</i>	16	64
Fosfomycin w/G6P	<i>Serratia</i> genus	32	256
Garenoxacin	<i>Stenotrophomonas maltophilia</i>	1	
Kanamycin	<i>Proteus mirabilis</i>	4	
Meropenem	<i>Proteus mirabilis</i>	0.5	
Ticarcillin	<i>Proteus mirabilis</i>	4	

Antimicrobial	Organism	Min MIC (µg/mL)	Max MIC (µg/mL)
Gram Positive			
Amikacin	<i>Staphylococcus epidermidis</i>	8	
	<i>Staphylococcus</i> other*	8	
Amoxicillin	<i>Staphylococcus aureus</i>	2	
	<i>Staphylococcus epidermidis</i>	2	
Ampicillin	<i>Staphylococcus aureus</i>		1
	<i>Enterococcus</i> other**	1	
Ampicillin-Sulbactam (f) (SXA)	<i>Staphylococcus epidermidis</i>	1	
Ampicillin-Sulbactam (SAM)	<i>Enterococcus</i> other**	2	
Cefotetan	<i>Staphylococcus epidermidis</i>	8	
Ceftriaxone	<i>Staphylococcus aureus</i>	2	
Cephalexin	<i>Staphylococcus epidermidis</i>	1	
Daptomycin	<i>Enterococcus casseliflavus</i>		
	<i>gallinarum</i>	2	
Lincomycin	<i>Enterococcus</i> other**	1	8
Netilmicin	<i>Staphylococcus epidermidis</i>	2	
Minocycline**	<i>Enterococcus</i> other**	1	32
	<i>Staphylococcus</i> other**	1	32
Norfloxacin	<i>Enterococcus</i> other**	2	
Penicillin	<i>Staphylococcus aureus</i>		1
	<i>Staphylococcus epidermidis</i>		1
	<i>Enterococcus</i> other**	2	
Piperacillin	<i>Staphylococcus epidermidis</i>	2	
Teicoplanin	<i>Staphylococcus epidermidis</i>	2	
Ticarcillin	<i>Staphylococcus epidermidis</i>	4	
	<i>Staphylococcus</i> other*	4	
Ticarcillin/Clavulanate	<i>Staphylococcus aureus</i>	4/2	
Trovaflaxacin	<i>Enterococcus faecium</i>	1	
Streptococci			
Daptomycin	Viridans group	0.25	

* Coagulase negative *Staphylococcus* other than *S. epidermidis*

** *Enterococcus* sp. other than *E. faecalis* and *E. faecium*

7.5 List of Reagents and Principles Employed in the BD Phoenix System

7.5.1 Gram Negative

SUBSTRATE NAME	CODE	PRINCIPLE
L-PHENYLALANINE-AMC	A_LPHET	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
4MU-N-ACETYL-BD-GLUCOS-AMINIDE	M_NAG	
L-GLUTAMIC ACID-AMC	A_LGTA	
L-TRYPTOPHAN-AMC	A_LTRY	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
L-PROLINE-AMC	A_LPROB	
L-ARGININE-AMC	A_LARGH	
ARGININE-ARGININE-AMC	A_ARARR	
GLYCINE-AMC	A_GLYB	
L-LEUCINE-AMC	A_LLEUH	
LYSINE-ALANINE-AMC	A_LYALD	
GLUTARYL-GLYCINE-ARGININE-AMC	A_GUGAH	
GLYCINE-PROLINE-AMC	A_GLPRB	
COLISTIN	C_CLST	Resistance to the antimicrobial agent results in a reduction of the resazurin based indicator.
POLYMYXIN B	C_PXB	
D-MANNITOL	C_DMNT	Utilization of a carbon source results in a reduction of the resazurin based indicator.
CITRATE	C_CIT	
ACETATE	C_ACT	
ADONITOL	C_ADO	
MALONATE	C_MLO	
ALPHA-KETOGLUTARIC ACID	C_KGA	
TIGLIC ACID	C_TIG	
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
L-PROLINE-NA	N_LPROT	Enzymatic hydrolysis of the colorless amide substrate releases yellow <i>p</i> -nitroaniline.
GAMMA-L-GLUTAMYL-NA	N_LGGH	
BIS (PNP) PHOSPHATE	P_BPHO	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow <i>p</i> -nitrophenol.
PNP-BD-GLUCOSIDE	P_BDGLU	

BETA-ALLOSE	R BALL	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
N-ACETYL-GALACTOSAMINE	R_NGA	
N-ACETYL-GLUCOSAMINE	R_NGU	
SORBITOL	R_DSBT	
SUCROSE	R_DSUC	
GALACTURONIC ACID	R_GRA	
MALTULOSE	R_MTU	
L-RHAMNOSE	R_LRHA	
BETA-GENTIOBIOSE	R_BGEN	
DEXTROSE	R_DEX	
D-GALACTOSE	R_DGAL	
D-FRUCTOSE	R_DFRU	
D-GLUCONIC ACID	R_DGUA	
D-MELIBIOSE	R_DMLB	
L-ARABINOSE	R_LARA	
METHYL-B-GLUCOSIDE	R_MBGU	
ORNITHINE	S_ORN	Utilization of ornithine results in pH rise and change in fluorescent indicator.
UREA	S_URE	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.

7.5.2 Gram Positive

SUBSTRATE NAME	CODE	PRINCIPLE
4MU-BD-CELLOBIOSIDE	M_BDCEL	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
L-ALANINE-AMC	A_LALT	
4MU-BD-GLUCOSIDE	M_BDGLU	
L-PROLINE-AMC	A_LPROB	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
L-PHENYLALANINE-AMC	A_LPHET	
L-TRYPTOPHAN-AMC	A_LTRY	
4MU-PHOSPHATE	M_PHOS	
METHIONINE-AMC	A_META	
4MU-AD-GLUCOSIDE	P_ADGLU	
ARGININE-ARGININE-AMC	A_ARARR	
GLYCINE-PROLINE-AMC	A_GLPRB	
4MU-BD-GLUCURONIDE	M_BDGLC	
L-LEUCINE-AMC	A_LLEUH	
4MU-N-ACETYL-BD-GLUCOS-AMINIDE	M_NAG	
L-ARGININE-AMC	A_LARGH	
4MU-PHOSPHATE (with Trehalose)	M_PHOT	
L-HISTIDINE-AMC	A_LHIST	
L-ISOLEUCINE-AMC	A_LISO	
4MU-BD-GALACTOSIDE	M_BDGAL	
COLISTIN	C_CLST	Resistance to the antimicrobial agent results in a reduction of the resazurin based indicator.
POLYMYXIN B	C_PXB	
D-GLUCONIC ACID	C_DGUA	Utilization of a carbon source results in a reduction of the resazurin based indicator.
3-METHYL GLUTARIC ACID	C_3MGA	
D-FRUCTOSE	C_DFRU	
IMINODIACETIC ACID	C_IMN	
ALPHA-KETOGLUTARIC ACID	C_KGA	
D-MANNITOL	C_DMNT	
3-METHYLADIPIC ACID	C_MAA	
THYMIDINE	C_THY	
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
ALANINE-ALANINE-PNA	N_ALALH	Enzymatic hydrolysis of the colorless amide substrate releases yellow <i>p</i> -nitroaniline.
L-PROLINE-PNA	N_LPROT	
VALINE-ALANINE-PNA	N_VAALA	
PNP-AD-GLUCOSIDE	P_PAGLU	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow <i>p</i> -nitrophenol.
PNP-PHOSPHATE	P_PHOL	

BETA-GENTIOBIOSE	R_BGEN	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
D-SUCROSE	R_DSUC	
MALTOTRIOSE	R_MTT	
N-ACETYL-GLUCOSAMINE	R_NGU	
D-TREHALOSE	R_DTRE	
D-TAGATOSE	R_DTAG	
MALTOSE	R_MAL	
DEXTROSE	R_DEX	
METHYL- α -D-GLUCOSIDE	R_MGP	
UREA	S_URE	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.
NITROCEFEN	L_NCF	Enzymatic hydrolysis of the β -Lactam ring results in a color change.

7.5.3 Streptococci Panel

SUBSTRATE NAME	CODE	PRINCIPLE
AMYGDALIN	R_AMY	Utilization of carbohydrate results in lower pH and change in indicator (Phenol red).
D-GALACTOSE	R_DGAL	
D-MANNITOL	R_DMTL	
D-RAFFINOSE	R_DRAF	
D-SORBITOL	R_DSBT	
D-TREHALOSE	R_DTRE	
DEXTRIN	R_DNX	
N-ACETYL-GLUCOSAMINE	R_NGU	
PHENYL GLUCOSIDE	R_PHG	
SALICIN	R_SAL	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow <i>p</i> -nitrophenol.
ONP-BD-GLUCOSIDE	O_BOGLU	
PNP-AD-GALACTOSIDE	P_ADGAL	
PNP-BD-CELLOBIOSIDE	P_CELB	
PNP-BD-GALACTOSIDE	P_GALB	
PNP-AD-GLUCOSIDE	P_PAGLU	
PNP-PHOSPHATE	P_PHOL	Enzymatic hydrolysis of the colorless amide substrate releases yellow <i>p</i> -nitroanilide.
ALANINE-ALANINE-PNA	N_ALALH	
VALINE-ALANINE-PNA	N_VAALA	
L-LYSINE-PNA	N_LLYSB	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
FLUORESCENT POSITIVE CONTROL	FLR_CTL	Utilization of a carbon source resulting in a reduction of the indicator (Resazurin based).
THYMIDINE	THY	
PULLULAN	PUL	
D-TREHALOSE	TRL	
D-LACTOSE	DLAC	

SUBSTRATE NAME	CODE	PRINCIPLE
LYSINE-AMC	A_LYSA	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
SERINE-TYROSINE-AMC	A_SETY	
L-CITRULLINE-AMC	A_LCTU	
L-PYROGLUTAMIC ACID-AMC	A_LPYR	
ISOLEUCINE-AMC	A_LISO	
L-TRYPTOPHAN-AMC	A_LTRY	
L-VALINE-AMC	A_LVAL	
ARGININE-ARGININE-AMC	A_ARARR	
LYSINE-ALANINE-AMC	A_LYALD	
ASPARAGINE-AMC	A_APGT	
L-ARGININE-AMC	A_LARGH	
L-HISTIDINE-AMC	A_LHIST	
ALANINE-AFC	Z_ALFT	
4MU-BD-CELLOBIOSIDE	M_BDCEL	
4MU-BD-GLUCOSIDE	M_BDGLU	
4MU-PHOSPHATE	M_PHOS	
4MU-AD-GLUCOSIDE	M_ADGLU	
4MU-BD-GLUCURONIDE	M_BDGLC	
4MU-N-ACETYL-BD-GLUCOSAMINE	M_NAG	
4MU-PHOSPHATE (with trehalose)	M_PHOT	
4MU-BD-GALACTOSIDE	M_BDGAL	
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.

7.5.4 Yeast Panel

SUBSTRATE NAME	CODE	PRINCIPLE
PNP-BD-GLUCOSIDE	P_BDGLU	Enzymatic hydrolysis of the colorless aryl substituted glycoside releases yellow <i>p</i> -nitrophenol.
PNP-AD-GLUCOSIDE	P_PAGLU	
ONP-BD-GLUCOSIDE	O_BOGLU	
L-SORBOSE	C_LSBO	Utilization of a carbon source results in a reduction of the resazurin based indicator.
DEXTROSE	C_DEX	
D-MANNITOL	C_DMNT	
D-SUCROSE	C_DSUC	
METHYL-AD-GLUCOPYRANOSIDE	C_MGP	
N-ACETYL-BD-GLUCOSAMINIDE	C_NAG	
DEXTROSE	R_DEX	Utilization of carbohydrate results in lower pH and change in indicator (phenol red).
D-FRUCTOSE	R_DFRU	
D-GALACTOSE	R_DGAL	
SUCROSE	R_DSUC	
D-TREHALOSE	R_DTRE	
MALTOTRIOSE	R_MTT	
ESCULIN	T_ESC	Hydrolysis of esculin results in a black precipitate in the presence of ferric ion.
FLUORESCENT NEGATIVE CONTROL	Z_FTST	Control to check for fluorescent interference.

FLUORESCENT POSITIVE CONTROL	FLR_CTL	Control to standardize fluorescent substrate results.
FLUORESCENT POSITIVE CONTROL	FLR_CTL	
GAMMA-L-GLUTAMYL-NA	N_LGGH	Enzymatic hydrolysis of the colorless amide substrate releases yellow <i>p</i> -nitroaniline.
L-PROLINE-PNA	N_LPROT	
ASPARAGINE-AMC	A_APGT	
L-ARGININE-AMC	A_LARGH	
L-GLUTAMINE-AMC	A_LGLNB	
L-TYROSINE-AMC	A_LTYO	
L-HISTIDINE-AMC	A_LHIST	
ORNITHINE-AMC	A_ORN	
THREONINE-AMC	A_THR	
HYDROXY PROLINE-AMC	A_LHYP	
4MU-N-ACETYL-BD-GLUCOSAMINE	M_NAG	
4MU-AD-GLUCOSIDE	M_ADGLU	
4MU-PHOSPHATE	M_PHOS	
LYSINE-ALANINE-AMC	A_LYALD	
GLYCINE-ARGININE-AMC	A_GLARH	Enzymatic hydrolysis of the amide or glycosidic bond results in the release of a fluorescent coumarin or 4-methylumbelliferone derivative.
ALANINE-AFC	Z_ALFT	
GLYCINE-AMC	A_GLYB	
L-CITRULLINE-AMC	A_LCTU	
L-GLUTAMIC ACID-AMC	A_LGTA	
L-VALINE-AMC	A_LVAL	
L-ALANINE-AMC	A_LALT	
L-PROLINE-AMC	A_LPROB	
L-TRYPTOPHAN-AMC	A_LTRY	
H-B-ALANINE-AMC	A_HBALT	
4MU-BD-CELLOBIOSIDE	M_BDCEL	
4MU-BD-GLUCOSIDE	M_BDGLU	
GLYCINE-PROLINE-AMC	A_GLPRB	
LYSINE-PROLINE-AMC	A_LYPRA	
BENZYL-L-CYSTEINE-AMC	A_BZLCY	
AMINO ACID	S_GTN	Utilization of the amino acid results in a change in fluorescence.
UREA	S_URE	Hydrolysis of urea and the resulting ammonia change results in pH rise and change in fluorescent indicator.

7.6 Taxa for ID/AST Determination

7.6.1 Gram Negative (0.5 McFarland)

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Achromobacter denitrificans</i>	AST	<i>Aeromonas trota</i>	AST
<i>Achromobacter piechaudii</i>	AST	<i>Aeromonas veronii</i> bv <i>sobria</i>	ID/AST
<i>Achromobacter</i> species	ID/AST	<i>Aeromonas veronii</i> bv <i>veronii</i>	ID/AST
<i>Achromobacter xylosoxidans</i>	AST	<i>Alcaligenes faecalis</i>	ID/AST
<i>Acinetobacter baumannii</i> ²	ID/AST	<i>Alcaligenes faecalis</i> ssp <i>faecalis</i>	AST
<i>Acinetobacter baumannii</i> / <i>calcoaceticus</i> complex ²	ID/AST	<i>Alcaligenes</i> species	AST
<i>Acinetobacter baumannii</i> / <i>haemolyticus</i>	AST	<i>Bergeyella zoohelcum</i>	ID
<i>Acinetobacter calcoaceticus</i>	AST	<i>Bordetella bronchiseptica</i>	ID
<i>Acinetobacter haemolyticus</i>	ID/AST	<i>Brevundimonas diminuta</i>	ID/AST
<i>Acinetobacter johnsonii</i>	AST	<i>Brevundimonas</i> species	AST
<i>Acinetobacter junii</i>	AST	<i>Brevundimonas vesicularis</i>	ID/AST
<i>Acinetobacter lwoffii</i>	ID/AST	<i>Burkholderia caryophylli</i>	AST
<i>Acinetobacter lwoffii</i> / <i>haemolyticus</i>	ID/AST	<i>Burkholderia cepacia</i> complex	ID/AST
<i>Acinetobacter radioresistens</i>	AST	<i>Burkholderia cepacia</i> / <i>Ralstonia pickettii</i>	ID/AST
<i>Acinetobacter</i> species	ID/AST	<i>Burkholderia gladioli</i>	ID/AST
<i>Actinobacillus lignieresii</i>	ID	<i>Burkholderia glathei</i>	AST
<i>Actinobacillus suis</i>	ID	<i>Burkholderia graminis</i>	AST
<i>Actinobacillus ureae</i>	ID	<i>Burkholderia multivorans</i>	AST
<i>Aeromonas allosaccharophila</i>	AST	<i>Burkholderia phenazinium</i>	AST
<i>Aeromonas caviae</i>	ID/AST	<i>Burkholderia pyrrocinia</i>	AST
<i>Aeromonas eucrenophila</i>	AST	<i>Burkholderia</i> species	AST
<i>Aeromonas hydrophila</i>	ID/AST	<i>Burkholderia</i> species/ <i>Ralstonia</i> species	ID/AST
<i>Aeromonas hydrophila</i> group	AST	<i>Cardiobacterium hominis</i>	ID
<i>Aeromonas jandaei</i>	AST	CDC group Vb-3	ID
<i>Aeromonas media</i>	AST	<i>Cedecea davisae</i>	ID/AST
<i>Aeromonas salmonicida</i>	AST	<i>Cedecea lapagei</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>achromogenes</i>	AST	<i>Cedecea neteri</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>masoucida</i>	ID/AST	<i>Cedecea</i> species	AST
<i>Aeromonas salmonicida</i> ssp <i>pectinolytica</i>	AST	<i>Cedecea</i> species 3	AST
<i>Aeromonas salmonicida</i> ssp <i>salmonicida</i>	ID/AST	<i>Cedecea</i> species 5	AST
<i>Aeromonas salmonicida</i> ssp <i>smithia</i>	ID/AST	<i>Chromobacterium violaceum</i>	ID
<i>Aeromonas schubertii</i>	ID/AST	<i>Chryseobacterium gleum</i>	ID/AST
<i>Aeromonas</i> species	AST	<i>Chryseobacterium indologenes</i>	ID/AST

* Taxa for AST Interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Chryseobacterium scopthalmum</i>	AST	<i>Escherichia fergusonii</i>	ID/AST
<i>Chryseobacterium species</i>	AST	<i>Escherichia hermannii</i>	ID/AST
<i>Citrobacter amalonaticus</i>	ID/AST	<i>Escherichia species</i>	AST
<i>Citrobacter braakii</i>	ID/AST	<i>Escherichia vulneris</i>	ID/AST
<i>Citrobacter farmeri</i>	ID/AST	<i>Ewingella americana</i>	ID
<i>Citrobacter freundii</i>	ID/AST	Gram-negative rod unidentified enteric	AST
<i>Citrobacter gillenii</i>	AST	Gram-negative rod unidentified nonfermenter	AST
<i>Citrobacter koseri</i>	ID/AST	<i>Grimontia hollisae</i>	ID
<i>Citrobacter murlinae</i>	AST	<i>Hafnia alvei</i>	ID/AST
<i>Citrobacter rodentium</i>	AST	<i>Hafnia alvei</i> group 1	AST
<i>Citrobacter sedlakii</i>	ID/AST	<i>Kingella denitrificans</i>	ID
<i>Citrobacter species</i>	AST	<i>Kingella kingae</i>	ID
<i>Citrobacter werkmanii</i>	ID/AST	<i>Klebsiella granulomatis</i>	AST
<i>Citrobacter youngae</i>	ID/AST	<i>Klebsiella oxytoca</i> ²	ID/AST
<i>Comamonas terrigena</i>	ID	<i>Klebsiella ozaenae</i>	ID/AST
<i>Comamonas testosteroni</i>	ID	<i>Klebsiella pneumoniae</i> ²	ID/AST
<i>Cosenzaea myxofaciens</i>	AST	<i>Klebsiella rhinoscleromatis</i>	ID/AST
<i>Cronobacter sakazakii</i> complex	ID/AST	<i>Klebsiella species</i>	AST
<i>Cupriavidus gilardii</i>	AST	<i>Kluyvera ascorbata</i>	ID/AST
<i>Cupriavidus pauculus</i>	ID/AST	<i>Kluyvera cryocrescens</i>	ID/AST
<i>Delftia acidovorans</i>	ID/AST	<i>Kluyvera georgiana</i>	AST
<i>Edwardsiella hoshinae</i>	ID/AST	<i>Kluyvera intermedia</i>	ID/AST
<i>Edwardsiella ictaluri</i>	ID/AST	<i>Kluyvera species</i>	AST
<i>Edwardsiella species</i>	AST	<i>Kosakonia cowanii</i>	AST
<i>Edwardsiella tarda</i>	ID/AST	<i>Leclercia adecarboxylata</i>	ID/AST
<i>Edwardsiella tarda</i> biogroup 1	AST	<i>Lelliottia amnigena</i>	AST
<i>Eikenella corrodens</i>	ID	<i>Lelliottia amnigena</i> biogroup 1	ID/AST
<i>Elizabethkingia meningoseptica</i>	ID/AST	<i>Lelliottia amnigena</i> biogroup 2	ID/AST
<i>Empedobacter brevis</i>	ID	<i>Leminorella grimontii</i>	ID
<i>Enterobacter aerogenes</i> ²	ID/AST	<i>Leminorella richardii</i>	ID
<i>Enterobacter asburiae</i>	ID/AST	<i>Mannheimia haemolytica</i>	ID
<i>Enterobacter cancerogenus</i>	ID/AST	<i>Methylobacterium extorquens</i>	ID
<i>Enterobacter cloacae</i> ²	ID/AST	<i>Moellerella wisconsensis</i>	ID/AST
<i>Enterobacter cloacae</i> ssp <i>dissolvens</i>	AST	<i>Moraxella (Branhamella) catarrhalis</i>	ID
<i>Enterobacter hormaechei</i>	ID/AST	<i>Moraxella species</i>	ID
<i>Enterobacter kobei</i>	AST	<i>Morganella morganii</i> ²	ID/AST
<i>Enterobacter nimipressuralis</i>	AST	<i>Morganella morganii</i> ssp <i>morganii</i>	AST
<i>Enterobacter species</i>	AST	<i>Morganella morganii</i> ssp <i>morganii</i> biogroup 1	AST
<i>Escherichia coli</i> ²	ID/AST	<i>Morganella morganii</i> ssp <i>sibonii</i>	AST
<i>Escherichia coli</i> serotype O111	AST	<i>Morganella species</i>	AST
<i>Escherichia coli</i> serotype O157	AST	<i>Myroides odoratimimus</i>	AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Myroides odoratus</i>	AST	<i>Pseudomonas luteola</i>	ID/AST
<i>Myroides odoratus/odoratimimus</i>	ID/AST	<i>Pseudomonas mendocina</i>	ID/AST
<i>Myroides species</i>	AST	<i>Pseudomonas monteilii</i>	AST
<i>Neisseria animaloris</i>	ID	<i>Pseudomonas oryzaehabitans</i>	ID/AST
<i>Neisseria zoodegmatidis</i>	ID	<i>Pseudomonas pertucinogena</i>	AST
<i>Ochrobactrum anthropi</i>	ID/AST	<i>Pseudomonas pseudoalcaligenes</i>	ID/AST
<i>Oligella ureolytica</i>	ID	<i>Pseudomonas pseudoalcaligenes</i> <i>ssp. pseudoalcaligenes</i>	AST
<i>Oligella urethralis</i>	ID	<i>Pseudomonas putida</i>	ID/AST
<i>Pantoea agglomerans</i>	ID/AST	<i>Pseudomonas species</i>	ID/AST
<i>Pantoea ananatis</i>	AST	<i>Pseudomonas stutzeri</i>	ID/AST
<i>Pantoea dispersa</i>	AST	<i>Pseudomonas veronii</i>	AST
<i>Pantoea species</i>	AST	<i>Rahnella aquatilis</i>	ID
<i>Pantoea stewartii</i>	AST	<i>Ralstonia pickettii</i>	ID/AST
<i>Pantoea stewartii ssp indologenes</i>	AST	<i>Ralstonia solanacearum</i>	AST
<i>Pantoea stewartii ssp stewartii</i>	AST	<i>Ralstonia species</i>	AST
<i>Paracoccus yeei</i>	ID	<i>Raoultella ornithinolytica</i>	ID/AST
<i>Pasteurella aerogenes</i>	ID	<i>Raoultella planticola</i>	AST
<i>Pasteurella multocida</i>	ID	<i>Raoultella species</i>	AST
<i>Pasteurella pneumotropica</i>	ID	<i>Raoultella terrigena</i>	AST
<i>Photobacterium damsela</i>	ID	<i>Rhizobium radiobacter</i>	ID
<i>Plesiomonas shigelloides</i>	ID	<i>Salmonella enterica ssp arizonae</i>	ID/AST
<i>Pluralibacter gergoviae</i>	ID/AST	<i>Salmonella enterica ssp diarizonae</i>	AST
<i>Pragia fontium</i>	ID	<i>Salmonella enterica ssp enterica</i> <i>serovar Choleraesuis</i>	ID/AST
<i>Proteus hauseri</i>	AST	<i>Salmonella enterica ssp enterica</i> <i>sv Gallinarum bv Gallinarum</i>	ID/AST
<i>Proteus mirabilis</i> ²	ID/AST	<i>Salmonella enterica ssp enterica</i> <i>sv Gallinarum bv Pullorum</i>	ID/AST
<i>Proteus penneri</i>	ID/AST	<i>Salmonella enterica ssp enterica</i> <i>sv Paratyphi A</i>	ID/AST
<i>Proteus species</i>	AST	<i>Salmonella enterica ssp enterica</i> <i>sv Typhi</i>	ID/AST
<i>Proteus vulgaris</i>	ID/AST	<i>Salmonella enterica ssp houtenae</i>	AST
<i>Proteus vulgaris/penneri</i>	ID/AST	<i>Salmonella enterica ssp indica</i>	AST
<i>Providencia alcalifaciens</i>	ID/AST	<i>Salmonella enterica ssp salamae</i>	AST
<i>Providencia heimbachae</i>	AST	<i>Salmonella species</i>	ID/AST
<i>Providencia rettgeri</i>	ID/AST	<i>Serratia entomophila</i>	AST
<i>Providencia rustigianii</i>	ID/AST	<i>Serratia ficaria</i>	ID/AST
<i>Providencia species</i>	AST	<i>Serratia fonticola</i>	ID/AST
<i>Providencia stuartii</i>	ID/AST	<i>Serratia grimesii</i>	AST
<i>Pseudomonas aeruginosa</i> ²	ID/AST	<i>Serratia liquefaciens</i>	ID/AST
<i>Pseudomonas alcaligenes</i>	AST	<i>Serratia marcescens</i> ²	ID/AST
<i>Pseudomonas fluorescens</i>	ID/AST	<i>Serratia odorifera</i>	AST
<i>Pseudomonas fluorescens/putida</i>	AST	<i>Serratia odorifera biogroup 1</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Serratia odorifera</i> biogroup 2	ID/AST	<i>Vibrio fluvialis</i>	ID
<i>Serratia plymuthica</i>	ID/AST	<i>Vibrio metschnikovii</i>	ID
<i>Serratia proteamaculans</i> ssp <i>proteamaculans</i>	AST	<i>Vibrio mimicus</i>	ID
<i>Serratia proteamaculans</i> ssp <i>quinovora</i>	AST	<i>Vibrio parahaemolyticus</i>	ID
<i>Serratia rubidaea</i>	ID/AST	<i>Vibrio vulnificus</i>	ID
<i>Serratia</i> species	AST	<i>Weeksella virosa</i>	ID
<i>Shewanella algae</i>	AST	<i>Yersinia aldovae</i>	AST
<i>Shewanella putrefaciens</i>	ID/AST	<i>Yersinia bercovieri</i>	AST
<i>Shewanella</i> species	AST	<i>Yersinia enterocolitica</i>	ID/AST
<i>Shigella boydii</i>	ID/AST	<i>Yersinia enterocolitica</i> group	AST
<i>Shigella dysenteriae</i>	ID/AST	<i>Yersinia frederiksenii</i>	ID/AST
<i>Shigella flexneri</i>	ID/AST	<i>Yersinia intermedia</i>	ID/AST
<i>Shigella sonnei</i>	ID/AST	<i>Yersinia kristensenii</i>	ID/AST
<i>Shigella</i> species	ID/AST	<i>Yersinia mollaretii</i>	AST
<i>Shimwellia blattae</i>	AST	<i>Yersinia pseudotuberculosis</i>	ID/AST
<i>Sphingobacterium multivorum</i>	ID/AST	<i>Yersinia rohdei</i>	AST
<i>Sphingobacterium multivorum</i> / <i>thalpophilum</i>	ID/AST	<i>Yersinia ruckeri</i>	ID/AST
<i>Sphingobacterium</i> species	AST	<i>Yersinia</i> species	AST
<i>Sphingobacterium spiritivorum</i>	ID/AST	<i>Yokenella regensburgei</i>	ID
<i>Sphingobacterium thalpophilum</i>	ID/AST		
<i>Sphingomonas paucimobilis</i>	ID/AST		
<i>Stenotrophomonas africana</i>	AST		
<i>Stenotrophomonas maltophilia</i> ²	ID/AST		
<i>Suttonella indologenes</i>	ID		
<i>Tatumella ptyseos</i>	ID		
<i>Vibrio alginolyticus</i>	ID		
<i>Vibrio cholerae</i>	ID		

* Taxa for AST Interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

7.6.2 Gram Negative (0.25 McFarland)

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Achromobacter</i> species	ID/AST	<i>Citrobacter youngae</i>	ID/AST
<i>Acinetobacter baumannii</i> / <i>calcoaceticus</i> complex	ID/AST	<i>Comamonas terrigena</i>	ID
<i>Acinetobacter haemolyticus</i>	ID/AST	<i>Comamonas testosteroni</i>	ID
<i>Acinetobacter Iwoffii</i>	ID/AST	<i>Cronobacter sakazakii</i> complex	ID/AST
<i>Actinobacillus lignieresii</i>	ID	<i>Cupriavidus pauculus</i>	ID/AST
<i>Actinobacillus suis</i>	ID	<i>Delftia acidovorans</i>	ID/AST
<i>Actinobacillus ureae</i>	ID	<i>Edwardsiella hoshinae</i>	ID/AST
<i>Aeromonas caviae</i>	ID/AST	<i>Edwardsiella ictaluri</i>	ID/AST
<i>Aeromonas hydrophila</i>	ID/AST	<i>Edwardsiella tarda</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>masoucida</i>	ID/AST	<i>Eikenella corrodens</i>	ID
<i>Aeromonas salmonicida</i> ssp <i>salmonicida</i>	ID/AST	<i>Elizabethkingia meningoseptica</i>	ID/AST
<i>Aeromonas salmonicida</i> ssp <i>smithia</i>	ID/AST	<i>Empedobacter brevis</i>	ID
<i>Aeromonas schubertii</i>	ID/AST	<i>Enterobacter aerogenes</i>	ID/AST
<i>Aeromonas veronii</i> bv <i>sobria</i>	ID/AST	<i>Enterobacter asburiae</i>	ID/AST
<i>Aeromonas veronii</i> bv <i>veronii</i>	ID/AST	<i>Enterobacter cancerogenus</i>	ID/AST
<i>Alcaligenes faecalis</i>	ID/AST	<i>Enterobacter cloacae</i>	ID/AST
<i>Bergeyella zoohelcum</i>	ID	<i>Enterobacter hormaechei</i>	ID/AST
<i>Bordetella bronchiseptica</i>	ID	<i>Escherichia coli</i>	ID/AST
<i>Brevundimonas diminuta</i>	ID/AST	<i>Escherichia fergusonii</i>	ID/AST
<i>Brevundimonas vesicularis</i>	ID/AST	<i>Escherichia hermannii</i>	ID/AST
<i>Burkholderia cepacia</i> complex	ID/AST	<i>Escherichia vulneris</i>	ID/AST
<i>Burkholderia gladioli</i>	ID/AST	<i>Ewingella americana</i>	ID
<i>Cardiobacterium hominis</i>	ID	<i>Grimontia hollisae</i>	ID
CDC group Vb-3	ID	<i>Hafnia alvei</i>	ID/AST
<i>Cedecea davisae</i>	ID/AST	<i>Klebsiella oxytoca</i>	ID/AST
<i>Cedecea lapagei</i>	ID/AST	<i>Klebsiella ozaenae</i>	ID/AST
<i>Cedecea neteri</i>	ID/AST	<i>Klebsiella pneumoniae</i>	ID/AST
<i>Chromobacterium violaceum</i>	ID	<i>Klebsiella rhinoscleromatis</i>	ID/AST
<i>Chryseobacterium gleum</i>	ID/AST	<i>Kluyvera ascorbata</i>	ID/AST
<i>Chryseobacterium indologenes</i>	ID/AST	<i>Kluyvera cryocrescens</i>	ID/AST
<i>Citrobacter amalonaticus</i>	ID/AST	<i>Kluyvera intermedia</i>	ID/AST
<i>Citrobacter braakii</i>	ID/AST	<i>Leclercia adecarboxylata</i>	ID/AST
<i>Citrobacter farmeri</i>	ID/AST	<i>Lelliottia amnigena</i> biogroup 1	ID/AST
<i>Citrobacter freundii</i>	ID/AST	<i>Lelliottia amnigena</i> biogroup 2	ID/AST
<i>Citrobacter koseri</i>	ID/AST	<i>Leminorella grimontii</i>	ID
<i>Citrobacter sedlakii</i>	ID/AST	<i>Leminorella richardii</i>	ID
<i>Citrobacter werkmanii</i>	ID/AST	<i>Mannheimia haemolytica</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST*
<i>Moellerella wisconsensis</i>	ID/AST	<i>Ralstonia pickettii</i>	ID
<i>Morganella morganii</i>	ID/AST	<i>Raoultella ornithinolytica</i>	ID/AST
<i>Myroides odoratus/odoratimimus</i>	ID/AST	<i>Rhizobium radiobacter</i>	ID
<i>Neisseria animaloris</i>	ID	<i>Salmonella enterica</i> ssp <i>arizonae</i>	ID/AST
<i>Neisseria zoodegmatidis</i>	ID	<i>Salmonella enterica</i> ssp <i>enterica</i> serovar <i>Choleraesuis</i>	ID/AST
<i>Ochrobactrum anthropi</i>	ID/AST	<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Gallinarum</i> bv <i>Gallinarum</i>	ID/AST
<i>Oligella ureolytica</i>	ID	<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Gallinarum</i> bv <i>Pullorum</i>	ID/AST
<i>Oligella urethralis</i>	ID	<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Paratyphi</i> A	ID/AST
<i>Pantoea agglomerans</i>	ID/AST	<i>Salmonella enterica</i> ssp <i>enterica</i> sv <i>Typhi</i>	ID/AST
<i>Paracoccus yeei</i>	ID	<i>Salmonella</i> species	ID/AST
<i>Pasteurella aerogenes</i>	ID	<i>Serratia ficaria</i>	ID/AST
<i>Pasteurella multocida</i>	ID	<i>Serratia fonticola</i>	ID/AST
<i>Pasteurella pneumotropica</i>	ID	<i>Serratia liquefaciens</i>	ID/AST
<i>Photobacterium damsela</i>	ID	<i>Serratia marcescens</i>	ID/AST
<i>Plesiomonas shigelloides</i>	ID	<i>Serratia odorifera</i> biogroup 1	ID/AST
<i>Pluralibacter gergoviae</i>	ID/AST	<i>Serratia odorifera</i> biogroup 2	ID/AST
<i>Pragia fontium</i>	ID	<i>Serratia plymuthica</i>	ID/AST
<i>Proteus mirabilis</i>	ID/AST	<i>Serratia rubidaea</i>	ID/AST
<i>Proteus penneri</i>	ID/AST	<i>Shewanella putrefaciens</i>	ID/AST
<i>Proteus vulgaris</i>	ID/AST	<i>Shigella boydii</i>	ID/AST
<i>Proteus vulgaris/penneri</i>	ID/AST	<i>Shigella dysenteriae</i>	ID/AST
<i>Providencia alcalifaciens</i>	ID/AST	<i>Shigella flexneri</i>	ID/AST
<i>Providencia rettgeri</i>	ID/AST	<i>Shigella sonnei</i>	ID/AST
<i>Providencia rustigianii</i>	ID/AST	<i>Sphingobacterium multivorum</i>	ID/AST
<i>Providencia stuartii</i>	ID/AST	<i>Sphingobacterium spiritivorum</i>	ID/AST
<i>Pseudomonas aeruginosa</i>	ID/AST	<i>Sphingobacterium thalpophilum</i>	ID/AST
<i>Pseudomonas fluorescens</i>	ID/AST	<i>Sphingomonas paucimobilis</i>	ID/AST
<i>Pseudomonas luteola</i>	ID/AST	<i>Stenotrophomonas maltophilia</i>	ID/AST
<i>Pseudomonas mendocina</i>	ID/AST	<i>Suttonella indologenes</i>	ID
<i>Pseudomonas oryzihabitans</i>	ID/AST	<i>Tatumella ptyseos</i>	ID
<i>Pseudomonas putida</i>	ID/AST	<i>Vibrio alginolyticus</i>	ID
<i>Pseudomonas stutzeri</i>	ID/AST	<i>Vibrio cholerae</i>	ID
<i>Rahnella aquatilis</i>	ID	<i>Vibrio fluvialis</i>	ID
<i>Vibrio metschnikovii</i>	ID	<i>Yersinia intermedia</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM NEGATIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Vibrio mimicus</i>	ID	<i>Yersinia kristensenii</i>	ID/AST
<i>Vibrio parahaemolyticus</i>	ID	<i>Yersinia pseudotuberculosis</i>	ID/AST
<i>Vibrio vulnificus</i>	ID	<i>Yersinia ruckeri</i>	ID/AST
<i>Weeksella virosa</i>	ID	<i>Yokenella regensburgei</i>	ID
<i>Yersinia enterocolitica</i>	ID/AST		
<i>Yersinia frederiksenii</i>	ID/AST		

^{*} Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

² Organism encountered in clinical trials and ≥ 20 strains available for antimicrobial/organism performance analysis.

7.6.3 Gram Positive (0.5 McFarland)

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST [*]	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST [*]
<i>Aerococcus urinae</i>	ID	<i>Corynebacterium matruchotii</i>	ID
<i>Aerococcus viridans</i>	ID	<i>Corynebacterium minutissimum</i>	ID
<i>Alloiooccus otitis</i>	ID	<i>Corynebacterium propinquum</i>	ID
<i>Arcanobacterium haemolyticum</i>	ID	<i>Corynebacterium pseudodiphtheriticum</i>	ID
<i>Bacillus cereus</i>	ID	<i>Corynebacterium pseudotuberculosis</i>	ID
<i>Bacillus circulans</i>	ID	<i>Corynebacterium renale</i>	ID
<i>Bacillus coagulans</i>	ID	<i>Corynebacterium striatum</i>	ID
<i>Bacillus licheniformis</i>	ID	<i>Corynebacterium ulcerans</i>	ID
<i>Bacillus megaterium</i>	ID	<i>Corynebacterium urealyticum</i>	ID
<i>Bacillus pumilus</i>	ID	<i>Corynebacterium xerosis</i>	ID
<i>Bacillus subtilis</i>	ID	<i>Dermabacter hominis</i>	ID
<i>Bacillus thuringiensis</i>	ID	<i>Demacoccus nishinomiyaensis</i>	ID
<i>Brevibacillus brevis</i>	ID	<i>Enterococcus asini</i>	AST
<i>Brevibacterium species</i>	ID	<i>Enterococcus avium</i>	ID/AST
<i>Cellulomonas turbata</i>	ID	<i>Enterococcus casseliflavus</i>	ID/AST
<i>Cellulosimicrobium cellulans</i>	ID	<i>Enterococcus casseliflavus/gallinarum</i>	ID/AST
<i>Corynebacterium amycolatum</i>	ID	<i>Enterococcus cecorum</i>	AST
<i>Corynebacterium amycolatum/minutissimum</i>	ID	<i>Enterococcus columbae</i>	AST
<i>Corynebacterium amycolatum/striatum</i>	ID	<i>Enterococcus dispar</i>	AST
<i>Corynebacterium bovis</i>	ID	<i>Enterococcus durans</i>	ID/AST
<i>Corynebacterium diphtheriae</i>	ID	<i>Enterococcus durans/faecium</i>	AST
<i>Corynebacterium jeikeium</i>	ID	<i>Enterococcus faecalis</i>	ID/AST
<i>Corynebacterium kutscheri</i>	ID	<i>Enterococcus faecalis/faecium</i>	AST

^{*} Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Enterococcus faecium</i>	ID/AST	<i>Leuconostoc citreum</i>	ID
<i>Enterococcus flavescens</i>	AST	<i>Leuconostoc lactis</i>	ID
<i>Enterococcus gallinarum</i>	ID/AST	<i>Leuconostoc mesenteroides</i> ssp <i>cremoris</i>	ID
<i>Enterococcus gilvus</i>	AST	<i>Leuconostoc mesenteroides</i> ssp <i>mesenteroides</i>	ID
<i>Enterococcus haemoperoxidus</i>	AST	<i>Leuconostoc pseudomesenteroides</i>	ID
<i>Enterococcus hirae</i>	ID/AST	<i>Listeria grayi</i>	ID
<i>Enterococcus hirae/faecium</i>	AST	<i>Listeria innocua</i>	ID
<i>Enterococcus malodoratus</i>	AST	<i>Listeria ivanovii</i>	ID
<i>Enterococcus moraviensis</i>	AST	<i>Listeria monocytogenes</i>	ID
<i>Enterococcus mundtii</i>	AST	<i>Listeria monocytogenes/innocua</i>	ID
<i>Enterococcus pallens</i>	AST	<i>Listeria welshimeri</i>	ID
<i>Enterococcus pseudoavium</i>	AST	<i>Lysinibacillus sphaericus</i>	ID
<i>Enterococcus raffinosus</i>	ID/AST	<i>Macrococcus caseolyticus</i>	ID
<i>Enterococcus raffinosus/avium</i>	AST	<i>Micrococcus luteus</i>	ID
<i>Enterococcus ratti</i>	AST	<i>Micrococcus lylae</i>	ID
<i>Enterococcus saccharolyticus</i>	AST	<i>Paenibacillus alvei</i>	ID
<i>Enterococcus species</i>	AST	<i>Paenibacillus macerans</i>	ID
<i>Enterococcus sulfureus</i>	AST	<i>Pediococcus acidilactici</i>	ID
<i>Erysipelothrix rhusiopathiae</i>	ID	<i>Pediococcus damnosus</i>	ID
<i>Gardnerella vaginalis</i>	ID	<i>Pediococcus dextrinicus</i>	ID
<i>Gemella haemolysans</i>	ID	<i>Pediococcus parvulus</i>	ID
<i>Gemella morbillorum</i>	ID	<i>Pediococcus pentosaceus</i>	ID
<i>Globicatella sanguinis</i>	ID	<i>Rhodococcus equi</i>	ID
<i>Helcococcus kunzii</i>	ID	<i>Rothia dentocariosa</i>	ID
<i>Kocuria kristinae</i>	ID	<i>Rothia mucilaginosa</i>	ID
<i>Kocuria rosea</i>	ID	<i>Staphylococcus arlettae</i>	AST
<i>Kocuria varians</i>	ID	<i>Staphylococcus aureus</i>	ID/AST
<i>Kytococcus sedentarius</i>	ID	<i>Staphylococcus aureus</i> ssp <i>anaerobius</i>	AST
<i>Lactococcus garvieae</i>	ID	<i>Staphylococcus aureus</i> ssp <i>aureus</i>	AST
<i>Lactococcus lactis</i> ssp <i>cremoris</i>	ID	<i>Staphylococcus aricularis</i>	ID/AST
<i>Lactococcus lactis</i> ssp <i>hordniae</i>	ID	<i>Staphylococcus capitis</i>	ID/AST
<i>Lactococcus lactis</i> ssp <i>lactis</i>	ID	<i>Staphylococcus capitis</i> ssp <i>capitis</i>	ID/AST
<i>Lactococcus plantarum</i>	ID	<i>Staphylococcus capitis</i> ssp <i>ureolyticus</i>	ID/AST
<i>Lactococcus raffinolactis</i>	ID	<i>Staphylococcus caprae</i>	ID/AST
<i>Leifsonia aquatica</i>	ID	<i>Staphylococcus carnosus</i>	ID/AST
<i>Staphylococcus carnosus</i> ssp <i>carnosus</i>	AST	<i>Staphylococcus pasteurii</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Staphylococcus carnosus</i> ssp <i>utilis</i>	AST	<i>Staphylococcus pettenkoferi</i>	ID/AST
<i>Staphylococcus chromogenes</i>	ID/AST	<i>Staphylococcus piscifermentans</i>	AST
<i>Staphylococcus chromogenes/hyicus</i>	ID/AST	<i>Staphylococcus pulvereri</i>	AST
<i>Staphylococcus coagulase-negative</i>	AST	<i>Staphylococcus saccharolyticus</i>	AST
<i>Staphylococcus coagulase-positive</i>	AST	<i>Staphylococcus saprophyticus</i>	ID/AST
<i>Staphylococcus cohnii</i>	ID/AST	<i>Staphylococcus saprophyticus</i> ssp <i>bovis</i>	AST
<i>Staphylococcus cohnii</i> ssp <i>cohnii</i>	ID/AST	<i>Staphylococcus saprophyticus</i> ssp <i>saprophyticus</i>	AST
<i>Staphylococcus cohnii</i> ssp <i>urealyticum</i>	ID/AST	<i>Staphylococcus schleiferi</i>	ID/AST
<i>Staphylococcus condimentii</i>	AST	<i>Staphylococcus schleiferi</i> ssp <i>coagulans</i>	ID/AST
<i>Staphylococcus delphini</i>	AST	<i>Staphylococcus schleiferi</i> ssp <i>schleiferi</i>	ID/AST
<i>Staphylococcus epidermidis</i>	ID/AST	<i>Staphylococcus sciuri</i>	ID/AST
<i>Staphylococcus equorum</i>	ID/AST	<i>Staphylococcus sciuri</i> ssp <i>caraticus</i>	AST
<i>Staphylococcus felis</i>	ID/AST	<i>Staphylococcus sciuri</i> ssp <i>rodentium</i>	AST
<i>Staphylococcus fleurettii</i>	AST	<i>Staphylococcus sciuri</i> ssp <i>sciuri</i>	AST
<i>Staphylococcus gallinarum</i>	ID/AST	<i>Staphylococcus simulans</i>	ID/AST
<i>Staphylococcus haemolyticus</i>	ID/AST	<i>Staphylococcus species</i>	AST
<i>Staphylococcus haemolyticus/lugdunensis</i>	ID	<i>Staphylococcus succinus</i>	AST
<i>Staphylococcus hominis</i>	ID/AST	<i>Staphylococcus succinus</i> ssp <i>casei</i>	AST
<i>Staphylococcus hominis</i> ssp <i>hominis</i>	AST	<i>Staphylococcus succinus</i> ssp <i>succinus</i>	AST
<i>Staphylococcus hominis</i> ssp <i>novobiosepticus</i>	AST	<i>Staphylococcus vitulinus</i>	ID/AST
<i>Staphylococcus hyicus</i>	ID/AST	<i>Staphylococcus warneri</i>	ID/AST
<i>Staphylococcus intermedius</i>	ID/AST	<i>Staphylococcus warneri/pasteuri</i>	AST
<i>Staphylococcus kloosii</i>	ID/AST	<i>Staphylococcus xylosus</i>	ID/AST
<i>Staphylococcus lentus</i>	ID/AST	<i>Streptococcus acidominimus</i>	ID
<i>Staphylococcus lugdunensis</i>	ID/AST	<i>Streptococcus agalactiae</i> (Strep. group B)	ID
<i>Staphylococcus lutrae</i>	AST	<i>Streptococcus anginosus</i>	ID
<i>Staphylococcus muscae</i>	AST	<i>Streptococcus bovis</i> (Strep. group D)	ID
<i>Streptococcus bovis</i> I (Strep. group D)	ID	<i>Streptococcus mitis</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Streptococcus bovis</i> II (Strep. group D)	ID	<i>Streptococcus mitis/pneumoniae</i>	ID
<i>Streptococcus canis</i>	ID	<i>Streptococcus mutans</i>	ID
<i>Streptococcus constellatus</i>	ID	<i>Streptococcus oralis</i>	ID
<i>Streptococcus cristatus</i>	ID	<i>Streptococcus parasanguinis</i>	ID
<i>Streptococcus dysgalactiae</i> ssp <i>dysgalactiae</i>	ID	<i>Streptococcus pneumoniae</i>	ID
<i>Streptococcus dysgalactiae</i> ssp <i>equisimilis</i>	ID	<i>Streptococcus porcinus</i>	ID
<i>Streptococcus equi</i>	ID	<i>Streptococcus pyogenes</i> (Strep. group A)	ID
<i>Streptococcus equi</i> ssp <i>equi</i>	ID	<i>Streptococcus salivarius</i>	ID
<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>	ID	<i>Streptococcus sanguinis</i>	ID
<i>Streptococcus equinus</i>	ID	<i>Streptococcus sobrinus</i>	ID
<i>Streptococcus gordonii</i>	ID	<i>Streptococcus uberis</i>	ID
<i>Streptococcus</i> group C/G (large colony)	ID	<i>Streptococcus vestibularis</i>	ID
<i>Streptococcus intermedius</i>	ID	<i>Trueperella pyogenes</i>	ID

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations.

7.6.4 Gram Positive (0.25 McFarland)

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Aerococcus urinae</i>	ID	<i>Gemella morbillorum</i>	ID
<i>Aerococcus viridans</i>	ID	<i>Globicatella sanguinis</i>	ID
<i>Alloiococcus otitis</i>	ID	<i>Helcococcus kunzii</i>	ID
<i>Dermacoccus nishinomiyaensis</i>	ID	<i>Kocuria kristinae</i>	ID
<i>Enterococcus avium</i>	ID/AST	<i>Kocuria rosea</i>	ID
<i>Enterococcus casseliflavus</i>	ID/AST	<i>Kocuria varians</i>	ID
<i>Enterococcus durans</i>	ID/AST	<i>Kytococcus sedentarius</i>	ID
<i>Enterococcus faecalis</i>	ID/AST	<i>Lactococcus lactis</i> ssp <i>cremoris</i>	ID
<i>Enterococcus faecium</i>	ID/AST	<i>Lactococcus lactis</i> ssp <i>hordniae</i>	ID
<i>Enterococcus gallinarum</i>	ID/AST	<i>Lactococcus plantarum</i>	ID
<i>Enterococcus hirae</i>	ID/AST	<i>Leuconostoc citreum</i>	ID
<i>Enterococcus raffinosus</i>	ID/AST	<i>Leuconostoc lactis</i>	ID
<i>Gemella haemolysans</i>	ID		
<i>Leuconostoc mesenteroides</i> ssp <i>mesenteroides</i>	ID	<i>Staphylococcus saprophyticus</i>	ID/AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations

GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*	GRAM POSITIVE TAXA ¹	ID, AST, ID/AST*
<i>Listeria innocua</i>	ID	<i>Staphylococcus schleiferi</i> ssp <i>coagulans</i>	ID/AST
<i>Listeria monocytogenes</i>	ID	<i>Staphylococcus schleiferi</i> ssp <i>schleiferi</i>	ID/AST
<i>Micrococcus caseolyticus</i>	ID	<i>Staphylococcus sciuri</i>	ID/AST
<i>Micrococcus luteus</i>	ID	<i>Staphylococcus simulans</i>	ID/AST
<i>Micrococcus lylae</i>	ID	<i>Staphylococcus vitulinus</i>	ID/AST
<i>Pediococcus acidilactici</i>	ID	<i>Staphylococcus warneri</i>	ID/AST
<i>Pediococcus damnosus</i>	ID	<i>Staphylococcus xylosus</i>	ID/AST
<i>Pediococcus dextrinicus</i>	ID	<i>Streptococcus agalactiae</i> (Strep. group B)	ID
<i>Pediococcus parvulus</i>	ID	<i>Streptococcus anginosus</i>	ID
<i>Pediococcus pentosaceus</i>	ID	<i>Streptococcus bovis</i> (Strep. group D)	ID
<i>Rothia mucilaginosa</i>	ID	<i>Streptococcus bovis</i> I (Strep. group D)	ID
<i>Staphylococcus aureus</i>	ID/AST	<i>Streptococcus bovis</i> II (Strep. group D)	ID
<i>Staphylococcus auricularis</i>	ID/AST	<i>Streptococcus constellatus</i>	ID
<i>Staphylococcus capitis</i>	ID/AST	<i>Streptococcus cristatus</i>	ID
<i>Staphylococcus caprae</i>	ID/AST	<i>Streptococcus equi</i>	ID
<i>Staphylococcus carnosus</i>	ID/AST	<i>Streptococcus gordonii</i>	ID
<i>Staphylococcus chromogenes</i>	ID/AST	<i>Streptococcus</i> group C/G (large colony)	ID
<i>Staphylococcus cohnii</i> ssp <i>cohnii</i>	ID/AST	<i>Streptococcus intermedius</i>	ID
<i>Staphylococcus cohnii</i> ssp <i>urealyticum</i>	ID/AST	<i>Streptococcus mitis</i>	ID
<i>Staphylococcus epidermidis</i>	ID/AST	<i>Streptococcus mutans</i>	ID
<i>Staphylococcus equorum</i>	ID/AST	<i>Streptococcus oralis</i>	ID
<i>Staphylococcus felis</i>	ID/AST	<i>Streptococcus parasanguinis</i>	ID
<i>Staphylococcus gallinarum</i>	ID/AST	<i>Streptococcus porcinus</i>	ID
<i>Staphylococcus haemolyticus</i>	ID/AST	<i>Streptococcus pyogenes</i> (Strep. group A)	ID
<i>Staphylococcus hominis</i>	ID/AST	<i>Streptococcus salivarius</i>	ID
<i>Staphylococcus hyicus</i>	ID/AST	<i>Streptococcus sanguinis</i>	ID
<i>Staphylococcus intermedius</i>	ID/AST	<i>Streptococcus sobrinus</i>	ID
<i>Staphylococcus kloosii</i>	ID/AST	<i>Streptococcus uberis</i>	ID
<i>Staphylococcus lugdunensis</i>	ID/AST	<i>Streptococcus vestibularis</i>	ID
<i>Staphylococcus pasteurii</i>	ID/AST		

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations

7.6.5 Streptococci

STREPTOCOCCI TAXA ¹	ID, AST, ID/AST [*]	STREPTOCOCCI TAXA ¹	ID, AST, ID/AST [*]
<i>Streptococcus acidominimus</i>	ID/AST	<i>Streptococcus ferus</i>	AST
<i>Streptococcus agalactiae</i> (Strep. group B)	ID/AST	<i>Streptococcus gallolyticus</i> ssp <i>macedonicus</i>	AST
<i>Streptococcus alactolyticus</i>	AST	<i>Streptococcus gordonii</i>	ID/AST
<i>Streptococcus alpha-hemolytic</i>	AST	<i>Streptococcus</i> group A (small colony)	AST
<i>Streptococcus anginosus</i>	ID/AST	<i>Streptococcus</i> group A (<i>Strep. pyogenes</i>)	AST
<i>Streptococcus anginosus</i> (previously <i>milleri</i>) group	ID/AST	<i>Streptococcus</i> group B (<i>Strep. agalactiae</i>)	AST
<i>Streptococcus beta-hemolytic</i> ACG (large colony)	AST	<i>Streptococcus</i> group C (large colony)	AST
<i>Streptococcus bovis</i> (Strep. group D)	AST	<i>Streptococcus</i> group C (small colony)	AST
<i>Streptococcus bovis</i> I (Strep. group D)	ID/AST	<i>Streptococcus</i> group C/G (large colony)	ID/AST
<i>Streptococcus bovis</i> II (Strep. group D)	ID/AST	<i>Streptococcus</i> group C/G (small colony)	AST
<i>Streptococcus canis</i>	ID/AST	<i>Streptococcus</i> group CFG (small colony)	AST
<i>Streptococcus constellatus</i>	ID/AST	<i>Streptococcus</i> group D (non-enterococcus)	AST
<i>Streptococcus constellatus</i> ssp <i>constellatus</i>	AST	<i>Streptococcus</i> group E	AST
<i>Streptococcus constellatus</i> ssp <i>pharyngis</i>	AST	<i>Streptococcus</i> group F	AST
<i>Streptococcus criceti</i>	AST	<i>Streptococcus</i> group G (large colony)	AST
<i>Streptococcus cristatus</i>	ID/AST	<i>Streptococcus</i> group G (small colony)	AST
<i>Streptococcus downei</i>	AST	<i>Streptococcus</i> group L	AST
<i>Streptococcus dysgalactiae</i>	AST	<i>Streptococcus hyointestinalis</i>	AST
<i>Streptococcus dysgalactiae</i> ssp <i>dysgalactiae</i>	ID/AST	<i>Streptococcus infantarius</i>	AST
<i>Streptococcus dysgalactiae</i> ssp <i>equisimilis</i>	ID/AST	<i>Streptococcus infantarius</i> ssp <i>coli</i>	AST
<i>Streptococcus dysgalactiae/canis</i>	ID/AST	<i>Streptococcus infantarius</i> ssp <i>infantarius</i>	AST
<i>Streptococcus equi</i>	ID/AST	<i>Streptococcus infantis</i>	AST
<i>Streptococcus equi</i> ssp <i>equi</i>	ID/AST	<i>Streptococcus iniae</i>	AST
<i>Streptococcus equi</i> ssp <i>zooepidemicus</i>	ID/AST	<i>Streptococcus intermedius</i>	ID/AST
<i>Streptococcus equinus</i>	ID/AST	<i>Streptococcus milleri</i> group	AST

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations

STREPTOCOCCI TAXA ¹	ID, AST, ID/AST*	STREPTOCOCCI TAXA ¹	ID, AST, ID/AST*
<i>Streptococcus mitis</i>	ID/AST	<i>Streptococcus salivarius</i>	ID/AST
<i>Streptococcus mitis</i> group	ID/AST	<i>Streptococcus salivarius</i> group	AST
<i>Streptococcus mitis/oralis</i>	AST	<i>Streptococcus salivarius</i> ssp <i>thermophilus</i>	AST
<i>Streptococcus mitis/pneumoniae</i>	ID/AST	<i>Streptococcus sanguinis</i>	ID/AST
<i>Streptococcus mutans</i>	ID/AST	<i>Streptococcus sanguinis</i> group	AST
<i>Streptococcus mutans</i> group	AST	<i>Streptococcus sobrinus</i>	ID/AST
<i>Streptococcus oralis</i>	ID/AST	<i>Streptococcus species</i>	AST
<i>Streptococcus parasanguinis</i>	ID/AST	<i>Streptococcus suis</i>	AST
<i>Streptococcus peroris</i>	AST	<i>Streptococcus uberis</i>	ID/AST
<i>Streptococcus pneumoniae</i>	ID/AST	<i>Streptococcus vestibularis</i>	ID/AST
<i>Streptococcus porcinus</i>	AST	<i>Streptococcus viridans</i> beta-hemolytic (small colony)	AST
<i>Streptococcus pyogenes</i> (Strep. group A)	ID/AST	<i>Streptococcus viridans</i> group	AST
<i>Streptococcus rattii</i>	AST		

* Taxa for AST interpretation may vary depending on the user-selected Interpretation Rule Set (Section 3.10.6 Custom Interpretation Rule Set Sub-Tab).

¹ Not all species encountered during clinical performance evaluations

7.6.6 Yeast

YEAST TAXA2	SABDX SABEM SABHI	CHOC COLSB TSASB
<i>Candida albicans</i>	√	√
<i>Candida apicola</i>	√	√
<i>Candida boidinii</i>		√
<i>Candida bracarensis</i>	√	
<i>Candida catenulata</i>	√	√
<i>Candida ciferrii</i>	√	√
<i>Candida dubliniensis</i>	√	√
<i>Candida firmetaria</i>	√	√
<i>Candida freyschussii</i>	√	√
<i>Candida glabrata</i>	√	√
<i>Candida guilliermondii</i>	√	√
<i>Candida guilliermondii</i> var <i>membranaefaciens</i>	√	√
<i>Candida haemulonii</i>	√	√
<i>Candida inconspicua</i>	√	√
<i>Candida kefyr</i>	√	√
<i>Candida krusei</i>	√	√
<i>Candida lipolytica</i>	√	√
<i>Candida lusitanae</i>	√	√

YEAST TAXA2	SABDX SABEM SABHI	CHOC COLSB TSASB
<i>Candida magnoliae</i>	√	√
<i>Candida melibiosica</i>	√	√
<i>Candida membranaefaciens</i>	√	√
<i>Candida norvegensis</i>	√	√
<i>Candida parapsilosis complex</i>	√	√
<i>Candida pararugosa</i>	√	√
<i>Candida pelliculosa</i>	√	√
<i>Candida pulcherrima</i>	√	√
<i>Candida rugosa</i>	√	√
<i>Candida sake</i>	√	√
<i>Candida sphaerica</i>	√	
<i>Candida tropicalis</i>	√	√
<i>Candida utilis</i>	√	√
<i>Candida viswanathii</i>	√	√
<i>Candida zeylanoides</i>	√	√
<i>Cryptococcus albidus</i>	√	√
<i>Cryptococcus humicola</i>	√	√
<i>Cryptococcus laurentii</i>		√
<i>Cryptococcus luteolus</i>	√	√
<i>Cryptococcus neoformans</i>	√	√
<i>Cryptococcus terreus</i>	√	
<i>Cryptococcus uniguttulatus</i>	√	√
<i>Exophiala dermatitidis</i>	√	√
<i>Exophiala species</i>	√	√
<i>Geotrichum species</i>	√	√
<i>Hortaea werneckii</i>	√	√
<i>Hyphopichia burtonii</i>	√	√
<i>Kloeckera species</i>	√	
<i>Magnusiomyces capitatus</i>	√	√
<i>Malassezia furfur complex</i>	√	
<i>Malassezia pachydermatis</i>	√	
<i>Malassezia sympodialis</i>	√	
<i>Millerozyma farinosa</i>	√	√
<i>Prototheca wickerhamii</i>	√	√
<i>Prototheca zopfii</i>	√	√
<i>Rhodotorula glutinis</i>	√	√
<i>Rhodotorula minuta</i>	√	√

YEAST TAXA2	SABDX SABEM SABHI	CHOC COLSB TSASB
<i>Rhodotorula mucilaginosa</i> var <i>mucilaginosa</i>	√	√
<i>Saccharomyces cerevisiae</i>	√	√
<i>Sporobolomyces salmonicolor</i>	√	√
<i>Trichosporon asahii</i>	√	√
<i>Trichosporon inkin</i>	√	
<i>Trichosporon loubieri</i>	√	
<i>Trichosporon mucoides</i>	√	√
<i>Trichosporon ovoides</i>		√
<i>Zygosaccharomyces bailii</i>	√	

¹ With an ID result of *Cryptococcus neoformans* and a substrate profile that may be indicative of *C. gattii*, the following message will print: The instrument produced an ID result of *Cryptococcus neoformans* with a profile that may be indicative of *C. gattii*. Recommend testing to rule out *C. gattii*.

² Not all species encountered during clinical performance evaluations.

8 – Glossaries

8.1 Glossary

AST	Antimicrobial Susceptibility Test(ing)
ATCC	American Type Culture Collection
BP	Blood plate
breakpoint	An interpretation of panel MIC data that which produces Susceptible, Intermediate, or Resistant result classifications. Breakpoints in the MIC data are established by the Clinical and Laboratory Standards Institute (CLSI) and other groups.
caddy	Accessory device used to transport inoculated panels to the instrument for loading.
carousel	The rotating drum which holds the 50 BD Phoenix panels and positions them for test readings, barcode readings, and loading or unloading.
CLSI	Clinical and Laboratory Standards Institute
CNA	Colistin Nalidixic Acid agar
demographic data	Accession information for a panel record
DTG	Drug Test Group
EMB	Eosin Methylene Blue agar
error station	Station that has sustained an optical or electro-mechanical error and has been blocked.
EUCAST	European Committee on Antimicrobial Susceptibility Testing
HE	Hektoen Enteric media agar
ID	Microbial Identification
ID/AST combination (combo) panel	The disposable device that contains all reagents needed for both ID and AST.
inoculation station	The inoculation station holds three BD Phoenix panels at the appropriate angle for optimal fill. The station also holds six broth tubes total, two per organism tested. One tube is for dilution of colony growth for Identification, the other for AST.

instrument test cycle	A complete test of all sample panels located in the carousel, resulting in color and/or fluorescence data values being recorded for each pertinent well of each panel.
MIC	Minimum Inhibitory Concentration; the lowest concentration of an antimicrobial which prohibits continued growth of the tested organism.
Needs Attention	Panels in the instrument's database that have encountered a condition that requires operator attention.
normalizer	A reference panel for use in the BD Phoenix M50 instrument. The Normalizer panel contains a matrix of visible light absorber and fluorescent material in panel-well format, which is used to correct individual well signals for losses occurring in the optical system.
orphan	A panel with a valid sequence number, but no associated accession number and isolate number.
panel carrier	The plastic carousel insert, which clips each BD Phoenix sample panel into place. The carousel contains 52 inserts, 26 in each of two tiers.
panel dataset	Each panel's set of color and fluorescence measurements, the panel's position identifier, test time stamps, and error flags are recorded for each test cycle throughout the test protocol. The test parameters determining each well's results are keyed by the individual barcode label signifying the panel type, and hence the type of test.
PEA	Phenylethyl Alcohol agar
position	The station. The physical location of the BD Phoenix panel within the instrument. This identifier includes instrument number, carousel tier letter, and numeric position on that tier.
panel presence detection threshold	Each inventory scan or test cycle, the instrument looks for both the barcode label of a panel in each carousel position and panel well data, color or fluorescence. If either is detected the instrument declares a panel logically present in that location. A panel without a valid barcode to provide panel type information will not be processed, but will be flagged as a Needs Attention candidate.
BD Phoenix AST indicator	An oxidation-reduction indicator used to signify microbial metabolism in the BD Phoenix panels. The indicator changes from blue to pink as initial reduction occurs. Further reduction causes the indicator to change from pink to colorless.
"rapid" results	AST Result obtained within 16 hours of panel inoculation.
related panel	Panels with the same accession number and isolate number are related.
Resistance Marker	Condition that is triggered when specific results indicate antimicrobial resistance. The action of some BDXpert rules is to trigger Resistance Markers; other rules may be called as a result of a specific Resistance Marker being triggered.
RGB	Red, Green and Blue. A shorthand representation of the visible light sources / wavelength regions used to interrogate the BD Phoenix panel.

sample	A specimen contained in a BD Phoenix panel. In practice, this would be a processed and resuspended dilution of microbiological growth from primary isolation culture in either ID diluent or AST broth which is then poured into the test panel.
SFM	Société Française de Microbiologie
SIR	Susceptible, Intermediate, or Resistant; refers to breakpoint AST categories. See also Breakpoint.
station	The instrument carousel is divided vertically into two tiers (A & B or C & D), each of which holds 26 panels. With one location occupied by a Normalizer panel, 25 locations per tier can accommodate test panels. This means that test panels can populate 50 total locations. Each location is assigned a tier letter and a number to determine the location on the tier. Indicator LEDs located behind each panel indicate station status (see Section 3 – Controls and Indicators).
test sequence number	A count of the number of readings taken by the instrument for a given sample well in a particular test panel, initialized to zero at the time of panel entry.
TSA	BD Trypticase Soy Agar
XLD	Xylose Lysine Deoxycholate agar

8.2 Supplemental ID Test Abbreviations

Test	Description	Test	Description
10C	Growth at 10 degrees Celsius	CAT	Catalase
42C	Growth at 42 degrees Celsius	CEL	Cellobiose
45C	Growth at 45 degrees Celsius	CIT	Citrate
50C	Growth at 50 degrees Celsius	COA	Coagulase
ACE	MBM + Acetate	DNA	Dnase
ALC	a-Lactose	ESC	Esculin
ANR	Anaerobic growth	FRU	Fructose
ARA	Arabinose	GAS	Gas from Glucose
ARG	Arginine	GEL	Gelatin
BE	Bile esculin	GLC	Beta glucuronidase
BSO	Bile solubility	GLU	Glucose
CAM	CAMP with <i>Staphylococcus aureus</i>	H ² S	Hydrogen Sulfide

Test	Description	Test	Description
HGN	Hemolysis - Gram Negative	OPS	Optochin susceptibility
HGP	Hemolysis - Gram Positive	ORN	Ornithine
HIP	Hippurate	OX	Oxidase
IND	Indole	PXR	Polymyxin Resistance
KCN	Growth in KCN	PXS	Polymyxin Susceptibility
LAC	Lactose	PYR	Pyrrolidonyl arylamidase
LYS	Lysine	RAF	Raffinose
MAC	Growth in MacConkey	SBT	Sorbitol
MAL	Maltose	SLT	Growth in 6.5% NaCl
MEL	Melibiose	SOR	Sorbose
MNS	Mannose	SUC	Sucrose
MNT	Mannitol	TRE	Trehalose
MOR	Morphology	URE	Urea
MOT	Motility	VAN	Vancomycin
MR	Methyl Red	VP	Voges Proskauer
NIT	Nitrate	XYL	Xylose
NVR	Novobiocin Resistance	YEL	Yellow/orange pigment

8.3 Symbol Glossary

Symbol	Meaning
	Manufacturer
	Catalog Number
	Authorized Representative in the European Community
	In Vitro Diagnostic Medical Device

9 – Replacement Parts

The following items may be ordered by contacting your local BD representative (See 11 – International Contacts).

Catalog Number	Item
443809	Barcode Scanner (external)
443866	System Software
441107	BD Phoenix™ Update Data (PUD) Software
443390	Barcode Scanner Stand
443842	Air Filter
443575	BD Phoenix™ M50 System User's Manual
443894	Quick Reference Guide
443431	Laser Printer
448984	Temperature Standard Panel

10 – Software Update Log

Whenever you receive a software update, please take a moment to log it below. This can assist you and BD personnel in identifying software revision levels, potential software problems, etc.

Date Received	Software Version	Date Installed	Installed By	Notes

11 – International Contacts

BD 7 Loveton Circle Sparks, MD 21152 USA Voice: (410) 316.4000 • Fax: (410) 527.0244 Technical Service and Support: 1.800.638.8663 Customer Service: 1.800.675.0908 www.bd.com	2100 Derry Road West Suite 100 Mississauga, ON Canada L5N 0B3 Voice: 905.288.6042 Fax: 905.288.6006
Monte Pelvoux 111 • 9th Floor Col. Lomas de Chapultepec 11000 Mexico D.F. Voice: 52 5 237 1200 • Fax: 52 5 237 1287	11 Rue Aristide Bergès BP4 38800 Le Pont de Claix France Voice: 33 476 68 36 36 • Fax: 33 476 68 34 95
Akasaka DS Building 5-26 Akasaka 8-chome Minato-ku, Tokyo 107-0052 Japan Voice: (81) 3 54138181 • Fax (81) 3 54138144	30 Tuas Avenue 2 Singapore 639461 Voice: (65) 8610633 • Fax: (65) 8601590
Rua Alexandre Dumas 1976 04717-004 Sao Paulo, S.P. Brazil Voice: (55) 11 5459833 • Fax: (55) 11 2478644	Australian Sponsor: Becton Dickinson Pty Ltd. 4 Research Park Drive Macquarie University Research Park North Ryde, NSW 2113 Australia

ATTACHMENT № 1

to OPERATIONAL DOCUMENTATION for Medical Device:

BD Phoenix Automated Microbiology System, variations: BD Phoenix 100,
BD Phoenix M50, with Accessories

Manufacturer:

Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152, USA

1. Product description

BD Phoenix Automated Microbiology System, variations: BD Phoenix 100, BD Phoenix M50, with Accessories

Variations:

I. BD Phoenix 100 Automated Microbiology System with Accessories.

Composition:

1. Base block
2. External barcode scanner
3. BD Phoenix 100 Quick Reference Guide (BD Phoenix - Quick Reference Guide, 448022)
4. BD Phoenix 100 software (PX SYSTEM SOFTWARE, 448035)

Accessories:

1. BD Phoenix starter kit (BD Phoenix - Starter Kit, 448040).
 - Inoculation Stations– 5 items;
 - Panel Caddies– 2 items;
 - Pipette tips - 960 items;
 - Thermometer – 1 item;
 - Spare air filters - 4 items;
 - External bar code scanner – 1 item;
 - Cable – 3 items;
 - Quick Reference Guide.
2. Table for BD Phoenix automated microbiology system (PX/PX AP TABLE, 441624).
3. BD Phoenix Inoculation Station (BD Phoenix - Inoculation Station, 448030).
4. BD Phoenix Panel Caddy (BD Phoenix - Panel Caddy, 448020)
5. BD Phoenix pipette stand (BD Phoenix - Pipette Stand, 448028).
6. BD Phoenix panel closures (BD Phoenix - Panel Closure, 448031) – not more than 200 items.
7. Thermometer (THERMOMETER PANEL, 448984).
8. Label panel carrier (LABEL PANEL CARRIER, 447111) - 1 sheet.
9. BD Phoenix pipette tips (BD Phoenix - Pipette Tips, 448037) – not more than 8 packages 204 tips each.

II. BD Phoenix M50 Automated Microbiology System with Accessories

Composition:

1. Base block
2. USB interconnect cable
3. Power cords for US, EU, China, UK – 4 pcs
4. Labels for instrument identification – 2 sets

Accessories:

1. PHOENIX M50 STARTER KIT, 443625:
 - All-in One PC: PC, UPS, cable, CD, bolts (4 items).
 - Instrument Power Cords (US/EU/China/UK) - 4 items
 - Monitor Mounting Arm Assembly & Instructions – 1 item
 - Socket Head Cap Screws – 2 pcs
 - Barcode Scanner
 - USB cables – 2 items
 - Barcode scanner stand
 - 25µL fixed volume pipette
 - 50µL fixed volume pipette
 - Pipette stand – 2 items

- 200µL pipette tips with filter – *1 pack of 8 boxes with 204 tips each*
- Phoenix Panel Transport Caddy – 2 items
- Phoenix Inoculation Station – 5 items
- Air Filter
- Blank USB Key

2. SOFTWARE PHX50, 443866

3. QUICK REFERENCE GUIDE PHX50 PKG, 443894

4. SCANNER BARCODE HANDHELD PACKAGED, 443809

5. MOUNTING ARM ERGOTRON PHX50, 443867

6. UPS SYNERGY SIII 2KVA ON-LINE, 257687

7. BATTERY PACK PAC SYN SIII 2KVA, 257688

8. Printer Brother HL54440D, 257613

9. KIT STACKING PHX50, 443910:

- bracket – 2 items,
- washer – 12 items,
- screw – 12 items,
- leveling foot – 2 items,
- drawing,
- anchor – 2 items,
- nuts hex – 2 items

10. SEISMIC KIT PHX50, 443918:

- bracket – 5 items,
- screw – 12 items,
- washer plain M5.3 X 10 X 1.0, SS – 12 items,
- bolt hex – 8 items,
- nut hex – 4 items,
- washer plain 656X 11.31 X.095 – 8 items,
- washer lock HSR. 625 – 8 items.

11. ADAPTER PLATE FX40 TO PHX50 SPARE, 443873:

- bracket – 6 items,
- screw – 12 items,
- washer – 8 items.

2. Information about the manufacturer of the medical product

MANUFACTURER

Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152, USA

MANUFACTURING SITES

Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152, USA.

3. Classification of the medical product

- Class of potential risk of medical product application: in compliance with the nomenclature specification: 2a
- GMDN code: 56747

4. Type of contact with the human body:

No contact with the body.

5. Purpose of the medical product

BD Phoenix Automated Microbiology System, variations: BD Phoenix 100, BD Phoenix M50, with

Accessories is an automated microbiology system, destined for quick identification (ID) and checking of antimicrobial susceptibility testing — AST of clinically significant bacteria. BD Phoenix automated microbiology system ensures quick results for the majority aerobic and some anaerobic gram-positive bacteria, as well as for the majority of aerobic and some anaerobic gram-negative bacteria, parasitizing on the human organism.

Indications for administration:

The product is not purposed for direct application on the patient.

The microbiology system is used only for *IN VITRO* microbiological diagnostics.

Contraindications:

No.

6. Sphere of application

Laboratory *IN VITRO* diagnostics and microbiology.

7. Description and technical specification of the product

System circuit diagram: See below (Fig. 1- 3)


Fig. 1 – Diagram of BD Phoenix 100 automated microbiology system incubation compartment


Fig. 2 – scheme of BD Phoenix 100 analyzer's incubator compartment

The main modules of BD Phoenix 100 Automated Microbiology System (hereinafter, Phoenix, Phoenix analyzer), which can be arbitrary derived, are:

1. Optical system and drive mechanism of the unit

Cylindrical turning gear and its motor are the only moving parts inside of the main incubator compartment. The turning gear unit is supported and rotated by three V-shaped bearing groups, which are coupled with V-shaped lips on the ring of turning gear. Usually turning gear rotates at a speed of 1.0 or 2.0 rpm. Speed, acceleration and deceleration of the motor are regulated by the CPU.

The turning gear is divided into four levels vertically, each of which operates as a relatively independent optical source and detection system for panels, placed at the level.

Each level has its own micro-controller to manage collection and transmission of data. Four microcontrollers communicates with the CPU via serial interface connection.

Visible radiation in red, green and blue spectral regions for each level is provided by the board of LED sources. Currents of the LED sources are programmed to compensate for the signal loss at the edges of panel due to displacement, and other factors. The source monitoring system averages the signal from two photodiodes of the board for tracking intensity of a source of visible light.

Fluorescent ultraviolet (UV) lamp that containing selected phosphor provides UV radiation. The active part of the lamp extends vertically through all four levels, providing continuous UV radiation accurate to 10% through the entire working length of the lamp. The second UV lamp is spare one.

2. Turning gear module

The turning gear is a frame structure made of aluminum rings and vertical ribs connected by screws forming a straight circular cylinder. The turning gear accommodates 104 panel holders on four levels. Each level has one normalizer panel and 25 sample panels.

Normally the turning gear rotates at a speed of 1.0 or 2.0 rpm depending on the operation. During the search for panel turning gear can rotate at speeds up to 10.0 rpm. The turning gear rotates counterclockwise. One rotation of the turning gear is used to perform pre-scanning to identify panels in the unit after closing the door. Complete test cycle takes seven minutes.

3 Incubation system

Incubation system maintains constant nominal temperature of the turning gear, panel holders, and the panels themselves at 35 °C. The system has a structure with forced air circulating convection and is equipped with a filter for dust removal. The system consists of a fan driven by a brushless DC motor, a ring-shaped heater with automatic shut-off in case of overheat, thermoformed system of input and output pipes, and a user-replaceable air filter made of polyester fibers.

4 Panel status and internal bar code scanner module

Panel status indication and its identification are performed on the same column of the unit directly below the turning gear, near the door. Panel status is displayed by tri-color LEDs mounted on the panel behind the turning gear, which glow through the fibers in the panel holders. With the door open one can see five columns of turning gear and five boards with four LEDs in each.

Four scanners reading bar codes of panels are arranged in the same column as the indicators. Each scanner reads the bar codes of the panels on the respective level of turning gear. Bar codes of panels are at the top of the base (bottom part) of the panel, when placing a panel in the unit they are facing inside.

5. Front panel and external bar code scanners

BD Phoenix 100 Automated Microbiology System is equipped by an embedded bar code scanner (bar code), which are located on the front side of the unit. (see fig 4).


Fig.4

For more convenient work BD Phoenix 100 Automated Microbiology System is also equipped by external bar code scanner plugging in through special port, which is located on the right side of the unit. Fig 5-6


Fig.5 – the right side of Phoenix 100


Fig.6 - External bar code scanner

These scanners can be used for reading inventory numbers of panels on the bar codes attached to them, and bar codes of own serial numbers of panels. Inventory numbers can be used to bind the sample identification information to specific panels in the unit.

6. Electronic control module

The electronic control module consists of the following two boards: IO interface, CPU (CPU) and memory, Ethernet, and liquid crystal display (LCD) driver. Other electronic components include user interface devices such as front panel bar code scanner, floppy disk drive, LCD display, special keypad, keyboard and loudspeaker for audio alarms.

Electronic control devices are primarily responsible for the following functions:

- planning and monitoring of panel testing;
- panel test data acquisition;
- turning on the indicators of panel (station) status;
- processing the data received from internal and external bar code scanners;
- interpretation of identification panels' data for identification of organisms;

- identification of antimicrobial drug susceptibility panels' data for receiving MIC and SIR results (in accordance with applicable rules);
- transmission of information to the LCD display;
- processing of keypad and keyboard commands;
- monitoring the state of door switch and lock;
- cabinet incubation control;
- interaction with external devices (such as EpiCenter workstation);
- alarm generation;
- performing built-in test (BIT).


Fig. 7 General view of BD Phoenix 100 analyzer


Fig. 8 BD Phoenix M50 Instrument Rear Fascia View


Fig. 9 BD Phoenix M50 Instrument Front Fascia View


Fig. 10 BD Phoenix M50 Instrument


Fig. 11 BD Phoenix M50 Instrument, shown in a stacked configuration


Fig. 12 BD Phoenix M50 Physical Block Diagram


Fig. 13 BD Phoenix M50 Functional Block Diagram


Fig. 14 BD Phoenix M50 Electrical Block Diagram


Fig. 15 BD Phoenix M50 Electronic Insulation Diagram

Table 1: Table of Electronic Insulation

Location	Insulation type	Working Voltage	Test Voltage	Creepage distance (mm)	Clearance (mm)
A	Double Insulation	240V~	2224V, 50Hz AC	6.0	3.0
B	Basic Insulation	240V~	1390V, 50Hz AC	3.0	1.5
C	Basic Insulation	240V~	1390V, 50Hz AC	3.0	1.5

The Phoenix M50 can be connected to the BD EpiCenter™ Plus Microbiology Data Management System. Connection of the Phoenix M50 instrument to Epicenter software has equivalent functionality to the connection between the Phoenix 100 instrument and Epicenter.

The EpiCenter system is a clinical data management system designed for use with certain BD instrument systems, including the Phoenix M50. The EpiCenter system provides advanced capabilities in data management, including extended demographic data, comprehensive reporting features including a powerful database query tool, bi-directional Laboratory Information System (LIS) communications and long-term data storage. Additionally, the EpiCenter system provides an Expert System that can manipulate and/or act upon AST results.

Table 2: Summary of Similarities and Differences

Feature/Functionality	Phoenix 100	Phoenix M50
Intended Use	The BD Phoenix™ Automated Microbiology System is intended for the rapid identification (ID) and antimicrobial susceptibility testing (AST) of clinically significant bacteria. The BD Phoenix system provides rapid results for most aerobic and facultative anaerobic gram-positive bacteria as well as most aerobic and facultative anaerobic gram-negative bacteria of human origin. The BD Phoenix system is also intended for the rapid identification of yeast and yeast-like organisms.	
Dimensions	105.9 cm H x 112.5 cm W x 74.9 cm D	535 mm H x 815 mm W x 765 mm D
Capacity	Up to 100 panels	Up to 50 panels
Ambient Operating Temperature	18 – 30°C	
Sample	Isolated colonies from pure culture.	Same
Inoculum preparation	Same inoculum preparation for both systems.	
Reagent system (Broths, panels, indicators)	Same reagent system for both systems	
Measurement system optics (Susceptibility results)	Same optical hardware for both systems	
Panel Transport Carousel Configuration	Four tiers	Two tiers
Panel Incubation Temperature	35°C ±1.5°C	
Algorithms	Same algorithms are present in both systems	
Performance	Equivalent performance in both systems	
Barcode Panel Tracking	Internal barcode scanner tracks the position of each panel in the instrument.	
Panel Entry	The door is opened and sample panels are entered. Light indicators illuminate to designate blocked stations. The user can put a panel in any non-blocked, empty station, which are not illuminated.	

8. Principle of work of the medical device

BD Phoenix 100 Automated Microbiology System can simultaneously perform up to 100 tests of identification and antimicrobial susceptibility test using Phoenix* combined panels. A maximum of 50 identification and antimicrobial susceptibility tests can be performed in the BD Phoenix M50 instrument at a time using BD Phoenix combination panels.

** Russian Federation FS registration certificate on the panel of BD Phoenix 100 analyzer: FSZ 2009/03983 of 16.03.2009*

Types of research: (ID) identification down to the genus and species of the organism. (AST) antimicrobial susceptibility test, method for producing best matching list of antibiotics to effectively treat the patient

The principles of analysis: Identification (ID) - colorimetric and fluorometric; Antibiotic susceptibility test (AST) - turbidity and color change indicator test.

Panel measuring and data processing take place every 20 minutes

4 sources of measurement are used

- Red spectrum
- Green spectrum
- Blue spectrum
- Ultraviolet spectrum

The only moving part in the unit is carrousel (which minimizes the possibility of failure). The carrousel has four independent rows of slots for panels with clinical specimens, with individual radiation sources

The unit utilizes two UV lamps: one of which is standby and switches on in case of the main lamp radiation quality deterioration.

Sealed molded polystyrene tray with automatic plating containing 136 micro-wells with dry reagents is consumable item for Phoenix system. Combined panel includes ID side with dry substrates for bacterial identification and AST side with different concentrations of antimicrobials, as well as standards of growth and fluorescence in corresponding wells. The Phoenix system uses an optimized colorimetric redox indicator for AST, as well as a number of colorimetric and fluorometric indicators for ID. AST broth is cationic (e.g., uses Ca ++ and Mg ++ ions) for optimization of susceptibility test performance.

Phoenix panel includes 51-well ID side and 85-well AST side. ID side contains 45 wells with dried biochemical substrates and 2 control fluorescent wells. AST side may hold up to 84 wells with dry antimicrobials and one growth control well. Panels only for ID, only for AST and combined ID/AST, are available. Unused wells are reserved for future use.

Phoenix panels are plated with standard inoculum. Suspensions of microorganisms should only be prepared using nephelometer BBL CrystalSpec, BD PhoenixSpec or BD Phoenix AP**. After plating the panels are placed in the unit and are constantly maintained at a temperature of 35 °C. The unit tests panels every 20 minutes: at the top of an hour, 20 minutes after the top of an hour and 40 minutes after the top of an hour - up to 16 hours if necessary. Information of Phoenix panels can only be read by a Phoenix device. Manual reading of Phoenix panel information is not allowed.

*** Russian Federation FS registration certificate on the BD PhoenixSpec nephelometer: FSZ: 2008/01073 of 18.02.2008*

Bacterial identification

Part of the Phoenix ID panel uses a range of traditional, chromogenic and fluorogenic biochemical tests to identify the organisms. To cover various types of activity in different taxa substrates based both on

growth and on the enzymes are used. The tests are based on the absorption and digestion of certain substrates by microorganisms, which are detected by different systems of indicators.

Acid formation is accompanied by the color change of phenol red indicator when the culture is capable of using hydrocarbon substrate. Chromogenic substrates are stained in yellow in the process of enzymatic hydrolysis of p-nitrophenyl or p-nitroanilide compounds. Enzymatic hydrolysis of fluorogenic substrates leads to the release of a fluorescent coumarin derivative. Organisms that use certain carbon sources recover an indicator based on resazurin. In addition, there are other tests to determine the ability of the organism to hydrolysis, degradation, recovery or other impact on substrate.

Antimicrobial susceptibility test

Phoenix AST method is a test with broth micro-dilution. The Phoenix system uses a redox indicator to determine the growth of organisms in the presence of an antimicrobial agent.

Continuous measurements of changes in the indicators, as well as the turbidity of bacteria are used to determine the growth of bacteria. Each AST panel configuration contains multiple antimicrobial agents with a broad range of dilution concentrations and one can be twice as high as the other. Identification of organisms is used in the interpretation of the values of MIC for each antimicrobial.

Required components for testing include: 1) Phoenix panels with lids 2) Phoenix ID broth; 3) Phoenix AST broth; 4) Phoenix AST indicator solution; 5) Phoenix AST-S broth (for use with panels Phoenix Strep only); 6) Phoenix AST-S indicator solution (for use with panels Phoenix Strep only); 7) Phoenix station; 8) Phoenix cassette for panels; 9) nephelometer BBL CrystalSpec, BD PhoenixSpec or BD Phoenix AP; 10) pipetting device with a capacity of 25 µl and sterile tips; 11) various laboratory gear.

Phoenix consumables are designed for use with the Phoenix device for identification of organisms (ID), antimicrobial susceptibility testing (AST) or both tests. The panel is positioned with the plating ports facing upward. Individual seed materials are added manually to the ID and AST ports. Seed materials penetrate deeply the panel through the winding canals, filling the panel wells as more liquid comes to the pad. Pad absorbs excessive seed material.

For manual filling the ports are equipped with caps. Air intake port is located near the panel cover separator to ensure adequate oxygen pressure in the panel during the test.

Principles of Phoenix AST tests for determining resistance markers

The following sections describe the principles of operation of Phoenix AST system for identifying resistance markers in Gram-negative and Gram-positive microorganisms, including the following: 1) The identification of ESBL forming in *Enterobacteriaceae*; 2) identification of *Vancomycin Resistant Enterococci* (VRE); 3) identification of high-level *aminoglycoside resistance in species of Enterococcus and Streptococcus* (HLAR); 4) identification of *methicillin resistance in staphylococci* (MRS); 5) Identification of *β-lactamase formation in the species of staphylococcus* (BL); 6) identification of *macrolides resistance in the species of Streptococcus* (MLSb and MEFF); 7) identification of high and low-level *penicillin resistance in S.pneumoniae* (HLPRSP) (LLPRSP).

Phoenix Extended Spectrum β-Lactamase (ESBL) test

Phoenix ESBL test is based on data from known ESBL antibiogram models. In order to develop Phoenix ESBL test there were selected strains of different microorganisms with known genotypes and phenotypes of β-lactamase from the *Enterobacteriaceae* families, including species *Escherichia coli*, *Klebsiella*, *Citrobacter*, *Enterobacter*, *Proteus*, and *Serratia*. Phoenix ESBL test is based on the principle of different reactions to the inhibitory effect of selected second or third generation cephalosporin in the presence or absence of a β-lactamase inhibitor, clavulanic acid.

Phoenix ESBL test principles are similar to the principles of a confirmatory test with CLSI ESBL.22 broth micro-dilution. Phoenix ESBL test is applied to the bacteria *E. coli*, *K. pneumonia*, and *K. oxytoca*. In addition, at the customer request, the test may be applied to other types of Enterobacteria if the formation of ESBL in it was documented in literature. In case of positive ESBL test result categorical interpretation of all penicillins, cephalosporins (except cephamycins) and aztreonam on the same Phoenix panel will be changed to stable.

Phoenix Vancomycin Resistant Enterococci (VRE)

Phoenix VRE test is based on the SIR interpretation of vancomycin. The breakpoint selected in the device configuration, is used for categorical interpretation. Phoenix VRE test was developed and optimized to be in accordance with the standard CLSI broth micro-dilution test.

Selecting a breakpoint different from CLSI can lead to suboptimal performance due to differences in the categorical interpretations. Only *Enterococcus faecalis* and *E. faecium* with acquired resistance (vanA or vanB) are considered positive results.

Phoenix High-Level Aminoglycoside Resistance (HLAR) test

Phoenix HLAR tests for *Enterococcus* species are based on the growth reaction in a single well with a high concentration of gentamicin or streptomycin. These tests were designed and optimized with respect to the CLSI broth micro-dilution tests and screening CLSI tests with agar.

Phoenix HLAR tests for *Streptococcus* species are based on the growth reaction in a single well containing gentamicin, kanamycin and streptomycin. These tests were designed and optimized using the recommended standard CLSI broth micro-dilution.

Phoenix Methicillin-Resistance in Staphylococci (MRS)

Phoenix MRS test is based on the interpretation of SIR oxacillin with Staphylococcus species. The breakpoint selected in the device configuration, is used for categorical interpretation. The test was designed and optimized with respect to the broth micro-dilution test under standard CLSI22, 23 and CLSI screening test with agar and oxacillin salt (BBL Oxacillin Screen Agar). Selection of breakpoint different from CLSI, can lead to suboptimal performance due to differences in categorical breakpoints in different standards. Using Phoenix MRS test allows testing all Staphylococcus specie. In case of positive MRS test results interpretations of all beta-lactam drugs in single Phoenix panel will be changed to stable.

Phoenix (β-lactamase (BL)) test

BL test, available in the Phoenix AST system, is a test of β- lactamase based on nitrocefin. The test based on nitrocefin, is a method of direct determination using ID side of Phoenix panel. The efficacy of this test was determined compared with the results of testing with Cefinase discs (BBL), used as a reference method. At present, these tests can only detect *Staphylococcus* species. In case of positive BL test results categorical interpretation of all labile penicillins on purified enzymes on the same Phoenix panels will be changed to stable.

Phoenix Macrolide Resistance in Streptococci (MLSb and MEFF) test

Phoenix Macrolide Resistance test is based on the interpretation of SIR erythromycin, clindamycin, and quinupristin/dalfopristin. The breakpoint selected in the device configuration, is used for categorical interpretation. Isolates of *Streptococcus* resistant to erythromycin and clindamycin and sensitive to quinupristin/dalfopristin, are phenotypically MLSb (macrolide/lincosamide/streptogramin B). For strains that are resistant to erythromycin and sensitive to clindamycin and quinupristin/dalfopristin the phenotype removing macrolides (Macrolide efflux - MEFF) or a MLSb inducing phenotype will be described. The formulation of clindamycin is tailored with respect to the goal of expanding definition of resistant phenotypes. This test was developed and optimized using the recommended standard CLSI broth micro-dilution.

Interpretation rules are defined according to Livermore and others. Phoenix Macrolide resistance test showed excellent correlation with genotypes *erm* and *mef* identified using molecular methods based on PCR.

Phoenix *mecA*-bound resistance marker for *Staphylococcus aureus* (*mecA*) species

Phoenix *mecA* test is used to predict *mecA*-bound resistance in *Staphylococcus aureus* species. The principle is similar to the principle of the disc diffusion test recommended by CLSI, which uses a cefoxitin disc (FOX) to predict *mecA*-bound resistance in *S. aureus* species. The efficiency of the test was determined by comparison with several PCR25 methods, as well as disc diffusion test. Using Phoenix *mecA* test the FOX breakpoints for *mecA*, used for defining resistance marker are set in the device. In case of detection of *mecA* resistance marker interpretations of all beta-lactam drugs in the single Phoenix panel is changed to stable and the Phoenix *mecA* resistance marker is set.

Phoenix Vancomycin Resistant *Staphylococcus aureus* (VRSA) test

Determination of Phoenix VRSA is based on SIR interpretation of vancomycin when testing *Staphylococcus aureus*. The breakpoint selected in the device configuration, is used for categorical interpretation. Phoenix VRSA test was developed and optimized to comply with the standard CLSI micro-dilution test and verified with well-known VRSA isolates. Selecting a breakpoint in different from

the one detected in CLSI M100-S15, may lower efficiency due to differences in the categorical interpretations. Only *Staphylococcus aureus* with true resistance (isolates containing the resistance mechanism, e.g. *vanA* gene) will be defined as VRSA. Strains of *S. aureus* with intermediate results in vancomycin (GISA/VISA) will be identified and described under different BDxpert rules. Gram-positive panel Phoenix AST determined resistance in VRSA *S. aureus* strains available at the time of comparative testing. The possibility of determining resistance of other *S. aureus* strains is unknown due to limited number of resistant strains that are available for comparative testing.

Phoenix High Level Penicillin Resistant *Streptococcus pneumonia* (HLPRSP) test

Phoenix Low Level Penicillin Resistant *Streptococcus pneumoniae* (LLPRSP) test

Phoenix HLPRSP test is based on the results of the minimum inhibitory concentration (MIC) of penicillin for *S. pneumoniae*. Penicillin MIC ≥ 4 indicates a high level of *S. pneumoniae* resistance to penicillin. Penicillin MIC in the range of from 0.125 to 2 indicates low level *S. pneumoniae* resistance to penicillin (LLPRSP).

9. System Performance

Definitions:

Essential Agreement (EA): essential agreement occurs when the results of MIC test of the Phoenix system and the test with micro-dilution of broth CLSI coincide or differ by no more than ± 1 dilution.

Categorical Agreement (CA): categorical agreement is observed when results of the MIC test of the Phoenix system and the test with micro-dilution of broth CLSI.

Efficiency of performance with Gram-negative bacteria

Identification of Gram-negative bacteria

The efficiency of identification of Gram-negative bacteria of the Phoenix system was measured using two internal studies. The values density of inoculum of 0.5 and 0.25 were tested with 721 strains (for the value 0.5) and 784 strains (for the value 0.25). Cultures of abdominal and non-abdominal type were evaluated by commercial and non-commercial methods.

The results of identification of Gram negative bacteria of the Phoenix system are given below:

	Density of inoculum (McFarland)	Coincidence	Lack of coincidence	Lack of identification
Level of genus	0.5	97.3%	2.1%	0.6%
	0.25	98.5%	1.0%	0.5%
Level of species	0.5	95.6%	3.6%	0.8%
	0.25	98.1%	1.4%	0.5%

Sensitivity of Gram-negative bacteria

Antimicrobial agent	Code	EA N	EA %	CA N	CA %
Amikacin	AN	974	94.3	974	99.2
Amoxicillin	AMX	725	95.9	839	96.9
Amoxicillin-Clavulanate	AMC	636	97.8	749	97.1
Amoxicillin-Clavulanate (f)	AXC	364	96.4	411	97.6
Ampicillin	AM	639	97.5	753	98.0
Ampicillin-Sulbactam	SAM	848	96.8	962	96.9
Ampicillin-Sulbactam (f)	SXA	1682	92.3	1682	94.5
Arbekacin	ARB	2083	94.1	2083	99.2
Aztreonam	ATM	1488	97.8	1431	98.5
Cefazolin	CZ	634	97.0	752	97.1
Cefdinir	CDR	633	92.7	747	94.5
Cefditoren	CDN	1270	97.2	1270	98.4

Antimicrobial agent	Code	EA N	EA %	CA N	CA %
Cefepime	FEP	1463	97.8	1463	98.6
Cefetamet	CAT	629	96.0	746	97.7
Cefixime	CFM	1454	96.7	1454	98.1
Cefmetazole	CMZ	608	94.7	718	96.4
Cefoperazone	CFP	854	95.1	972	97.7
Cefoperazone-Sulbactam	SCP	1921	93.1	1921	96.0
Cefotaxime	CTX	849	96.9	970	97.2
Cefotetan	CTT	627	96.2	748	96.4
Cefotiam	CFT	1428	94.9	1428	97.6
Cefoxitin	FOX	628	96.7	748	97.6
Cefpirome	CPO	846	96.2	964	97.0
Cefpodoxime	CPD	616	94.5	737	97.0
Cefsulodin	CFS	55	98.2	56	98.2
Ceftazidime	CAZ	848	95.9	969	97.5
Ceftibuten	CTB	590	90.5	708	95.9
Ceftizoxime	ZOX	854	97.3	971	97.2
Ceftriaxone	CRO	2416	96.1	2416	91.6
Cefuroxime	CXM	623	97.0	744	98.5
Cephalexin	CN	319	97.8	366	96.5
Cephalothin	CF	613	98.0	731	98.5
Chloramphenicol	C	978	97.7	978	98.7
Ciprofloxacin	CIP	977	99.4	977	99.7
Colistin	CL	467	97.9	467	98.7
Ertapenem*	ETP	1632	97.7	1265	99.4
Fosfomycin	FF	446	93.0	446	96.9
Garenoxacin	GRN	1977	98.1	1977	99.4
Gatifloxacin	GAT	752	99.5	752	99.6
Gemifloxacin	GEM	2096	98.6	2096	99.0
Gentamicin	GM	973	95.8	973	99.8
Imipenem	IPM	1545	94.8	1483	99.0
Isepamicin	ISP	468	93.8	468	99.4
Kanamycin	K	735	95.8	735	99.9
Levofloxacin	LVX	972	99.6	972	100.0
Lomefloxacin	LOM	976	99.1	976	99.6
Mecillinam	MEC	345	95.7	345	94.8
Meropenem	MEM	1450	97.0	1394	98.3
Minocycline	MI	2094	93.8	2094	98.3
Moxalactam	MOX	2063	96.9	2063	98.3
Moxifloxacin	MXF	746	98.5	747	99.7
Nalidixic Acid	NA	750	94.0	750	99.1
Netilmicin	NET	974	96.7	974	99.3
Nitrofurantoin	FM	744	98.4	744	98.8
Norfloxacin	NOR	976	98.9	976	99.6
Ofloxacin	OFX	971	99.3	971	99.7
Pefloxacin	PEF	469	98.5	469	99.4
Piperacillin	PIP	860	94.7	973	97.1
Piperacillin-Tazobactam	TZP	856	92.9	970	96.2
Temocillin	TEM	1410	96.1	1410	98.9
Tetracycline	TE	975	95.6	975	98.5
Ticarcillin	TIC	859	94.8	973	97.8
Ticarcillin-Clavulanate	TIM	534	92.7	589	94.6
Tigecycline	TGC	1428	97.1	1110	96.7
Tobramycin	NN	977	94.2	977	99.4
Trimethoprim	TMP	752	96.0	752	99.6
Trimeth-Sulfa (DIN)	STG	463	97.8	463	97.8
Trimethoprim-Sulfamethoxazole	SXT	976	95.5	976	97.3

*Ability to Detect Resistance Unknown

Efficiency of performance with Gram-positive bacteria

Identification of Gram-positive bacteria

The efficiency of identification of Gram-positive bacteria of the Phoenix system was measured using two internal studies. The values density of inoculum of 0.5 and 0.25 were tested with 696 strains (for the value 0.5) and 755 strains (for the value 0.25). The results were compared with the results of commercial and non-commercial methods.

Results of identification of Gram-positive bacteria of the Phoenix system are given below:

	Density of inoculum (McFarland)	Coincidence	Lack of coincidence	Lack of identification
Level of genus	0.5	99.1%	0.3%	0.6%
	0.25	99.6%	0.0%	0.4%
Level of species	0.5	95.4%	3.9%	0.7%
	0.25	98.0%	1.6%	0.4%

Sensitivity of Gram-positive bacteria

Antimicrobial agent	Code	EA N	EA %	CA N	CA %
Amikacin	AN	487	95.7	487	95.7
Amoxicillin	AMX	395	90.6	659	97.1
Amoxicillin-Clavulanate	AMC	180	97.8	446	97.5
Amoxicillin-Clavulanate (f)	AXC	397	96.7	664	97.6
Ampicillin	AM	402	94.0	667	98.7
Ampicillin-Sulbactam	SAM	179	96.7	445	97.1
Ampicillin-Sulbactam (f)	SXA	1449	94.6	1449	95.2
Arbekacin	ARB	973	96.2	973	99.9
Azithromycin	AZM	702	95.7	702	97.7
Cefaclor	CEC	182	91.8	449	95.6
Cefazolin	CZ	180	97.2	441	97.5
Cefdinir	CDR	182	94.5	434	96.8
Cefditoren	CDN	944	96.2	944	97.8
Cefepime	FEP	181	97.2	446	97.5
Cefmetazole	CMZ	174	96.6	430	97.2
Cefoperazone	CFP	184	95.7	447	97.8
Cefotaxime	CTX	185	96.2	446	97.8
Cefotetan	CTT	183	94.5	439	96.6
Cefotiam	CFT	965	91.5	965	94.0
Cefoxitin	FOX	184	96.2	445	97.5
Cefozopran	CFZ	1130	89.8	1130	89.6
Cefpirome	CPO	183	100.0	446	97.5
Cefpodoxime	CPD	185	91.4	451	95.6
Ceftaroline*	CPT	1313	94.5	866	99.7
Ceftazidime	CAZ	183	95.1	448	97.1
Ceftizoxime	ZOX	184	93.5	447	96.2
Ceftriaxone	CRO	184	95.1	447	97.3
Cefuroxime	CXM	184	94.6	447	98.2

Antimicrobial agent	Code	EA N	EA %	CA N %	CA %
Cephalexin	CN	185	90.3	452	97.4
Cephalothin	CF	182	97.8	441	97.3
Chloramphenicol	C	705	92.8	705	95.3
Ciprofloxacin	CIP	686	96.4	686	97.4
Clarithromycin	CLR	485	96.3	485	98.1
Clindamycin	CC	703	98.2	703	98.3
Daptomycin*	DAP	1361	97.6	1361	99.8
Doxycycline	D	1211	96.5	1211	99.8
Ertapenem	ETP	670	91.6	670	95.4
Erythromycin	E	474	95.4	474	95.8
Fosfomycin	FF	472	97.0	472	98.3
Fusidic Acid	FA	477	96.7	477	99.0
Garenoxacin	GRN	1212	97.1	1212	99.0
Gatifloxacin	GAT	481	98.8	481	100.0
Gemifloxacin	GEM	1400	96.2	1400	97.4
Gentamicin	GM	487	98.4	487	99.2
Gentamicin-Syn	GMS	NA	NA	198	97.5
Gentamicin-Syn (SFM)	GMF	NA	NA	87	100.0
Imipenem	IPM	185	100.0	448	97.5
Inducible Macrolide Resistance (iMLSB) Test	ECC	NA	NA	295	97.6
Kanamycin	K	463	92.2	464	97.4
Kanamycin Synergy	KS	501	100.0	501	96.6
Levofloxacin	LVX	698	96.9	698	98.6
Lincomycin	L	696	98.3	696	98.7
Linezolid	LZD	474	99.4	474	100.0
Lomefloxacin	LOM	486	97.7	486	99.2
Meropenem	MEM	185	98.4	447	97.5
Minocycline	MI	1448	98.5	1448	99.2
Moxalactam	MOX	935	95.9	935	96.0
Moxifloxacin	MXF	486	99.0	486	99.6
Mupirocin	MUP	938	96.6	938	99.7
Mupirocin High Level	MUH	968	100.0	968	99.9
Netilmicin	NET	488	93.9	488	96.3
Nitrofurantoin	FM	707	98.3	707	98.3
Norfloxacin	NOR	674	96.0	674	97.8
Ofloxacin	OFX	487	99.2	488	99.8
Oxacillin	OX	449	95.8	449	97.8
Pefloxacin	PEF	706	97.9	706	98.0
Penicillin	P	401	96.0	662	98.0
Piperacillin	PIP	382	81.7	639	98.0
Piperacillin-Tazobactam	TZP	386	100.0	649	98.3
Pristinamycin	PR	695	93.2	695	95.5
Quinupristin-Dalfopristin	SYN	704	97.2	704	97.6
Rifampin	RA	700	84.0	700	88.3
Streptomycin-Syn	STS	NA	NA	198	99.0
Streptomycin-Syn (SFM)	STF	NA	NA	87	98.8
Teicoplanin	TEC	671	97.6	671	99.4
Telithromycin	TEL	898	97.3	898	98.9
Tetracycline	TE	707	96.3	707	98.6
Ticarcillin	TIC	395	95.4	659	98.2
Ticarcillin-Clavulanate	TIM	396	97.5	656	97.9
Tigecycline*	TGC	1356	96.8	1044	100
Tobramycin	NN	473	97.0	473	97.3
Trimethoprim	TMP	479	95.4	479	97.9
Trimeth-Sulfa (DIN)	STG	645	91.0	646	92.7
Trimethoprim-Sulfamethoxazole	SXT	654	89.5	654	93.0
Vancomycin	VA	702	99.3	702	99.9

*Ability to Detect Resistance Unknown

Efficiency of performance with streptococci (with panels SMIC/ID and SMIC)

Identification of streptococci (only SMIC/ID)

The efficiency of identification of streptococci of the Phoenix system was measured using the internal study. The results for the 655 cultures were compared to the results of commercial and non-commercial methods.

The results of identification of streptococci of the Phoenix system are given below:

	Coincidence	Lack of coincidence	Lack of identification
Level of species	96.3%	2.4%	1.2%

The sensitivity of streptococci

The efficiency of the Phoenix system for streptococci AST was evaluated in four clinical laboratories using clinical cultures. The results of AST obtained from the panels of tests with micro-dilution of broth were compared, prepared in accordance with requirements of CLSI standard. Discrepant results were **arbitrated by duplicate repeat testing in both systems.**

Antimicrobial agent	Code	EA N	EA %	CA N	CA %
Amoxicillin	AMX	1576	99.2	1576	99.7
Amoxicillin-Clavulanate	AMC	1564	97.9	1564	99.3
Ampicillin	AM	1569	98.0	1569	99.4
Cefepime	FEP	1571	97.6	1571	99.7
Cefotaxime	CTX	1578	98.5	1578	99.6
Ceftriaxone	CRO	1579	98.5	1579	99.8
Cefuroxime	CXM	1581	97.7	1581	98.8
Chloramphenicol	C	1587	97.4	1587	99.6
Clindamycin	CC	1591	94.0	1591	97.4
Daptomycin	DAP	1566	93.1	1566	99.9
Ertapenem*	ETP	1585	97.9	1585	99.6
Erythromycin	E	1577	93.2	1577	98.2
Garenoxacin*	GRN	1515	95.5	1515	99.7
Gatifloxacin	GAT	1587	95.1	1587	99.5
Gemifloxacin*	GEM	1592	98.9	1592	99.7
Gentamicin-Syn	GMS	1032	99.6	1032	99.6
Imipenem	IPM	1581	98.0	1581	99.6
Kanamycin-Syn	KS	1032	99.6	1032	99.6
Levofloxacin	LVX	1595	97.4	1595	99.8
Linezolid	LZD	1586	97.1	1586	99.4
Meropenem	MEM	1579	97.4	1579	99.7
Moxifloxacin*	MXF	1590	95.8	1590	99.6
Ofloxacin	OFX	1594	98.2	1594	98.7
Penicillin	P	1585	93.9	1585	98.5
Pristinamycin*	PR	1583	94.8	1583	99.8
Quinupristin-dalfopristin*	SYN	1581	97.3	1581	99.7
Streptomycin-Syn	STS	1031	99.9	1031	99.9
Teicoplanin*	TEC	1593	99.8	1593	99.9
Telithromycin*	TEL	1574	95.9	1574	98.5
Tetracycline	TE	1589	93.3	1589	97.6
Trimethoprim-Sulfamethoxazole	SXT	1008	93.8	1008	98.3
Vancomycin*	VA	1588	98.2	1588	99.9

*Ability to Detect Resistance Unknown

11. Restrictions of the method

To precise special cases of restrictions for antimicrobials or microorganisms, see package insert, supplied together with the panel.

General characteristics. Gram stain test is necessary for selection of the respective types of Phoenix panels. Precise identification and/or AST cannot be performed without these tests.

Use only thoroughly selected bacterial colony from one of the recommended initial mediums. Use of mixed colonies could result in inaccurate identification and/or interpretation of AST.

If the density of the device inoculum (for the used type of panel) is 0.5, value of inoculum density 0.50-0.60 McFarland should be used. To measure inoculum density, it is possible to use BBL CrystalSpec or BD PhoenixSpec Nephelometer.

If inoculum density (for the used type of panel) is 0.25, value of inoculum density 0.20 – 0.30 McFarland should be used. To measure inoculum density of this range, only BD PhoenixSpec Nephelometer can be used.

Information from Phoenix panels is read only using Phoenix device. Visual interpretation of Phoenix panels is not allowed. Any attempts to manually interpret panel results may result in incorrect identification and/or inaccuracy in AST interpretation.

The panels should be placed into Phoenix instrument within 30 minutes after inoculation.

To obtain the most reliable results, it is recommended to reseed quality control microorganisms at least two times in two consecutive days into TSA with 5% Sheep Blood before being used in the Phoenix system.

Identification

Unique panel environment, joined with reduced incubation time, may entail Phoenix panel reaction, different from those obtained using conventional biochemical medium.

Testing for susceptibility to antimicrobials. Indicator solution AST may be added to test tube with AST broth and stored in dark place for up to 8 hours before use. Use only one drop of indicator solution AST from the bottle with a pipette. If more than one drop will be added by mistake, throw the test tube away and use a new test tube with AST broth. After adding AST indicator solution into the test tube with AST broth, mix them turning the tube upside down. DO NOT AGITATE. Agitation may result in appearance of bubbles in AST broth, which may entail incorrect inoculation of Phoenix panel during inoculation. Some microorganisms, included into Taxa identification, are not included into AST database due to low probability of their appearance. In case of revelation of such microorganisms, «Organism not included in the AST database, perform alternate method» message will appear.

For some combinations of microorganisms/antimicrobials, absence of stable strains hinders determination of any categories of results, besides «susceptible». For strains, causing results in «nonsusceptible» category, identification of microorganisms and the results of test of susceptibility to antimicrobials should be confirmed. After that, the cultures should be preserved and submitted for analysis to the authorized laboratory for confirmation of the result with the use of CLSL broth microdilution.

Use of CHROMagar Orientation may cause incorrect values of susceptibility when testing erythromycin with gram-positive microorganisms. The results of tests of antimicrobials susceptibility should be confirmed using Trypticase Soy Agar with 5% Sheep Blood.

This type of medium (“chocolate agar”) should not be used for identification of streptococcus with panels SMIC/ID. “Chocolate” agar may be used only for testing of streptococcus susceptibility.

Use of Columbia agar with 5% horse blood may result in significantly higher MIC value for SXT with kinds of Streptococcus, which may result in incorrect susceptibility value. The results of tests of susceptibility to antimicrobials should be confirmed using soya agar Trypticase with 5% Sheep Blood.

NOTE

The results of MIC for cultures without identification of species and kind in «Final ID» field (Final identification) can not be valid.

In case of modification of «Final ID» (Final identification) values of the cultures from that not revealed in list of AST Taxa to the declared MIC result and interpretation will be based on the final identification, which was declared in the list of AST Taxa.

12. System specification

BD Phoenix 100 Automated Microbiology System with Accessories:

Physical characteristics

Height 105.9 cm

Width 112.5 cm

Depth 74.9 cm

Distance to the wall (from the right, from the left and from behind) 10.2 cm

Distance to the wall (from front) 45.7 cm

Weight (without inoculation) 151.4 kg

Weight (full) 159.6 kg

The stand should bear the weight of 227.3 kg

Principles of the analysis:

ID identification: colorimetric and fluorometric

Determination of antibiotic sensitivity (AST): measurement of cloudiness and of color change indicator.

Requirements to electric mains

Input voltage: 110 – 117 V AC or 220 – 240 V DC $\pm 10\%$

Input current: 8A

Line input frequency: 47 - 53 Hz or 57 - 63 Hz

Capacity: 1.03 kW

Heat energy: 3 515 BTU/h

Environmental requirements

Storage conditions

Temperature: 17.8 – $> 65.0\text{ }^{\circ}\text{C}$

Humidity: Relative humidity: 10 – 90% without condensation

Exploitation conditions

Temperature: 17.8 – $> 30\text{ }^{\circ}\text{C}$

Humidity: Relative humidity: 20 – 90% without condensation

Location: Even surface, no direct sunlight and heat sources.

In the areas with elevated seismological activity, it is strictly recommended to use the appropriate fastening.

Installation category: II pollution rate: 2 in accordance with IEC 664.

Optical characteristics: wavelength peak values

From 428 to 623 nm, visible spectrum

367 nm, UV-excitation

Pass band: 410 – 641 nm, UV light

Presence of display:

LCD

BD Phoenix M50 Automated Microbiology System with Accessories:

Physical characteristics

Height 53.5 cm

Width 81.5 cm

Depth 76.5 cm

Clearance-Right 54.6 cm

Clearance-Left, Rear 7.6 cm

Clearance-Front 45.7 cm

Weight-Empty 53.5 kg

Weight-Full 57.6 kg

Principles of the analysis:

ID identification: colorimetric and fluorometric

Determination of antibiotic sensitivity (AST): measurement of cloudiness and of color change indicator.

Environmental requirements

Non-Operating Storage

Temperature: -17.8 – > 65.0 °C

Humidity: Relative humidity: 20 – 80% non-condensing

Operating conditions

Temperature: 18 – > 30 °C

Humidity: Relative humidity: 20 – 80% non-condensing

Locations: Level surface, no direct sunlight and no direct heat.

Altitude 2000 meters

Noise at 1 m ≤ 55 dBA using ANSI Type 2 sound meter

Electrical requirements

Input Voltage 100-240 VAC

Input Current 6 Amp

Input Line Frequency 50/60 Hz

Instrument Heat Output at Input Voltage =120 VAC

Ambient Temperature	BTU/hr
18 °C	648
24 °C	464
30 °C	287

Optical specifications: peak wavelengths

428 nm to 623 nm, visible spectrum

375 nm, UV Excitation

410 nm – 640 nm Bandpass, UV Emission

Presence of display:

LCD

13. Other functional specification of medical device

Range of optical density

- The ID inoculum starting point is 0.50 McFarland, or approximately 1.5×10^8 cfu/ml (E. coli)
- The low inoculum ID is 0.25 McFarland, or approximately 0.75×10^8 cfu/ml (E. coli)
- The AST inoculum is approximately 5×10^5 cfu/ml

Maximum permissible error

The visible spectrum optical transmission of the normalizer is 25%; the panel yields blue fluorescence in a band centered about 440 nm when illuminated with 367 nm UV light. The fluorescence signal is equivalent to that produced by 7.5 – 10.0 micromolar AMC at pH=8.0. The Normalizer serves as a reference panel for adjustment of the instrument's optical detection system. It is used for adjustment of the LED's (Red, Green and Blue) and fluorescence lamps.

The Normalizer Panel is always located in Position Zero (0) of each tier. Routine access to the carousel does not present the Normalizer Panel to the user when the carousel door is opened. The Normalizer Panel is constrained in its carrier to prevent inadvertent removal by the user. Access to the Normalizer Panel is controlled via LCD keypad selection of instrument diagnostics.

Wavelength of radiation dispersion of blue, green, red panels

The wavelength characteristics of the source LED's is presented in **Table 1-1.**

LED	Peak wavelength	Bandwidth at half- maximum
-----	-----------------	----------------------------

	(nm)	intensity (nm)
Red	623	15
Green	525	45
Blue	428	65

Source LED Characteristics Table 1-1

Time of full analyze

After a door-opening event, the instrument reads panel barcode labels and then performs a scan using the red LED's to determine if panels are present, and to locate or "map" well positions. Panel readings are made on the hour and at 20 and 40 minutes after the hour.

The instrument test sequence is listed below:

1. Lock the door
2. Wait for home flag to synchronize CPU processing with the tiers.
3. Set test cycle start time to the current time
4. Set the drum speed to 1 RPM.
5. Send start-scan command to the tiers
6. Delay 10 seconds while the tiers perform dark readings
7. Turn the UV bulb on
8. Wait 2 revolutions for UV Warm-up.
9. Take UV reading – One drum revolution
10. Turn UV bulb off
11. Set drum speed to 2 RPM
12. Turn on red LED's
13. Wait 1 revolution for Red Warm-up.
14. Take Red reading – One revolution
15. Turn off Red LED's
16. Turn on Green LED's
17. Wait 1 revolution for Green Warm-up.
18. Take Green readings – One revolution
19. Turn off Green LED's
20. Set drum speed to 1 RPM.
21. Turn on Blue LED's
22. Wait 1 revolution for Blue Warm-up.
23. Take Blue readings – One revolution
24. Turn off Blue LED's
25. Save this test cycle's start time as the start time of the last cycle
26. Unlock door
27. Post-test cycle processing

Productivity of analyzer

A maximum of 100 identification and antimicrobial susceptibility tests can be performed in the Phoenix instrument at a time using Phoenix combination panels.

The instrument tests panels every 20 minutes: on the hour, at 20 minutes past the hour, and again at 40 minutes past the hour for up to 16 hours if necessary. Each panel will take up to 16 hours to complete, however, the panel could be completed before the 16 hours. It depends on the organism and the time the ID and/or sensitivity takes to complete.

14. Description of parameters of BD Phoenix Software

Parameter	Description	Conformance BD Phoenix 100 Software to required parameter
1.1 . Function stability	Ability to ensure continuation of the programme after the deviations caused by hardware failures, errors in the input data and service errors	Conforms
1.2 . Working capacity	The program's ability to function within the specified modes and levels of information processed in accordance with policies in the absence of hardware failures	Conforms
2.1 . Structural properties	Organization of all interrelated parts of the program into a cohesive one using logical structures "sequence", "option", "repetition"	Conforms
2.2 .Simple design	Construction of the modular structure of the program the most rational way from the point of view of perception and understanding	Conforms
2.3 . Visibility	The existence and representation in most easily understanding way of source modules PS, full description of them in the relevant program documents	Conforms
2.4 . Repeatability	The use of standard designs or features within PS	Conforms
3.1 . Learnability	Presentation of the program papers and programs in the form of facilitating the understanding of the logic of operation of the program as a whole and its parts	Conforms
3.2 . Availability of operational program documents	Clarity, transparency and completeness of the descriptions of the user interaction with the program in operational program documents	Conforms
3.3 . Ease of use and service	Conformance the data processing and forms of result presentation to the nature of the tasks	Conforms
4.1 . Level of automation	The level of automation of data processing functions, taking into account the effectiveness of the	Conforms

	functional structure of the program in terms of the interaction with the user and the use of computing resources	
4.2 . Time efficiency	The program's ability to perform the specified actions in an interval of time that meets the specified requirements	Conforms
4.3 . Resource intensity	The minimum required computing resources and the number of service personnel for operation PS	Conforms
5.1 . Flexibility	The ability to use PS in a variety of applications	N/A
5.2 . Mobility	The ability to use PS without significant additional work on computers	N/A
5.3 . Modifiability	Simplicity of making the necessary adjustments and modifications to the program during operation	Conforms
6.1 . Realization fullness	Realization fullness of specified functions of the PS and the adequacy of their descriptions in the program documentation	Conforms
6.2 . Compatibility	Clear consistent description and use of identical objects, functions, terms, definitions, identities etc in different parts of the program documents and text programs	Conforms
6.3 . Logical correctness	Functional and program conformance the data processing to system requirements by tasks	Conforms
6.4 . Certainty	Completeness of checking the possible methods of program work during testing	Conforms

Table 3: Comparison of Software Features


Software Feature/Function	Phoenix 100	Phoenix 50
Operating System used on measurement system computer	<i>Phoenix 100 measurement subsystem computers use a real time operating system (VRTX RTOS) for high level control, the low level computers do not use an operating system.</i>	<i>The Phoenix M50 measurement subsystem computers use a real time operating system (smxBlaze RTOS)</i>
	Functional operation of the measurement system within the Phoenix M50 was accomplished through use of code from the Phoenix 100. Changes to the code were only made when the physical hardware environment required change. Leverage of the Phoenix 100 code assured that the operation of the Phoenix M50 measurement system was functionally equivalent.	
Operating system used on Instrument application computer	<i>Real time operating system (VRTX RTOS)</i>	<i>Microsoft Windows Embedded 7</i>
	Functional equivalence of the algorithmic processing performed with the M50 environment was accomplished by leveraging code from the Phoenix 100. Calculations of results is accomplished through Phoenix 100 code. The visualization of the information to the user on the Phoenix M50 is physically different, however, feature wise the information presented provides all the functionality found within the Phoenix 100.	
Expert System	BDXpert (Available on stand-alone Phoenix instrument, or on EpiCenter if using Phoenix/EpiCenter Configuration)	
User Interface	<i>User interface software on board instrument.</i>	<i>Primary user interface moved to accompanying PC. Instrument will allow users to load and unload panels and provide panel status updates.</i>

BD Phoenix automated microbiology system includes the following components:

- I. BD Phoenix 100 Automated Microbiology System with Accessories:
Composition:
 1. Base block


2. External barcode scanner

External bar code scanner supplying for more convenient work BD Phoenix 100 Automated Microbiology System is plunged in through special port, which is located on the right side of the unit. (see detail description above).

3. BD Phoenix 100 Quick Reference Guide (BD Phoenix - Quick Reference Guide, 448022)

4. BD Phoenix 100 software (PX SYSTEM SOFTWARE, 448035)

Special software for work with BD Phoenix 100 automated microbiology system.
On carrier of FDD 1.44 mb 3.5 type.

II. Accessories

1. BD Phoenix starter kit (BD Phoenix - Starter Kit, 448040).
The starter kit is purposed for automated microbiology system startup.

- Inoculation Stations– 5 items;
Plastic holders for comfortable work with panels.

- Panel Caddies– 2 items;
Caddies for transportation of panels.

- Pipette tips - 960 items;
Single-use polymeric pipette tips are used for samples preparation.


- Thermometer – 1 item;
Alcohol thermometer for placement in the analyzed incubation compartment. See detail description further. Is used for checking temp sensor readings.

- Spare air filters - 4 items;
HEPA - filters for cleaning of air, supplied to the incubation compartment.


- External bar code scanner – 1 item;
Additional scanner for panels plugging in through special port, which is located on the right side of the unit (see detail description above).

- Cable – 3 items;
Network power cable for the automated microbiology system.

- Quick Reference Guide.

2. Table for BD Phoenix automated microbiology system (PX/PX AP TABLE, 441624).
The automated microbiology system table is purposed for equipment installation.

3. BD Phoenix Inoculation Station (BD Phoenix - Inoculation Station, 448030).


Description: see above.

4. BD Phoenix Panel Caddy (BD Phoenix - Panel Caddy, 448020)

Description: see above.

5. BD Phoenix pipette stand (BD Phoenix - Pipette Stand, 448028).

Plastic stand for comfortable placement of pipettes on the working table.

6. BD Phoenix panel closures (BD Phoenix - Panel Closure, 448031) – not more than 200 items.


Package: 200 pcs.

Polymeric closures for panel closure.

Are used for panel sealing after preparation of clinical sample for analysis.

7. Thermometer (THERMOMETER PANEL, 448984).

Type thermometer: Alcohol thermometer

Range of measurement: from 30 to 37 °C.

This measurement range is stipulated by normal working temperature interval: 35°C ±1.5.

Graduating mark 0.2 °C


measurement inaccuracy: ±0.2 °C

Dimension (length*width*height, mm): 16*80*15

Weight: 35 g


Intended for use: for checking temp sensor readings. For temperature control in the incubation compartment of the analyzer.

8. Label panel carrier (LABEL PANEL CARRIER, 447111) - 1 sheet.

Adhesive stickers for carousel cells marking.

9. BD Phoenix pipette tips (BD Phoenix - Pipette Tips, 448037) – not more than 8 packages 204 tips each.


8 packages 204 tips each.

Description: see above.

II. BD Phoenix M50 Automated Microbiology System with Accessories:

Composition:

1. Base block

2. USB interconnect cable

3. Power cords for US, EU, China, UK – 4 items


4. Labels for instrument identification – 2 sets

To identify and number the instruments

Accessories:


1. PHOENIX M50 STARTER KIT, 443625:
 - All-in One PC: PC, UPS, cable, CD, bolts (4 items).

- Instrument Power Cords (US/EU/China/UK) - 4 items


- Monitor Mounting Arm Assembly & Instructions - 1 item
To connect the Instrument with All-in-one PC

- Socket Head Cap Screws - 2 items. To attach Mounting Arm to the instrument

- Barcode Scanner


- USB cables – 2 *items*.
To connect peripherals


- Barcode scanner stand
To place the barcode scanner.


- 25 μ L fixed volume pipette, 50 μ L fixed volume pipette & pipette stand (2 items)
To execute inoculation

- 200 μ L pipette tips with filter - *1 pack of 8 boxes with 204 tips each*

- Phoenix Panel Transport Caddy – 2 items
To transport panels to the instrument

- Phoenix Inoculation Station – 5 items

- Air Filter

- Blank USB Key


2. SOFTWARE PHX50, 443866
USB memory stick with system application software and user's guides


3. QUICK REFERENCE GUIDE PHX50 PKG, 443894
Hard copy user's guide book


4. SCANNER BARCODE HANDHELD PACKAGED, 443809

5. MOUTING ARM ERGOTRON PHX50, 443867

6. UPS SYNERGY SIII 2KVA ON-LINE, 257687
Protects key electronic equipment from power disruptions

7. BATTERY PACK PAC SYN SIII 2KVA, 257688
Used together with UPS SYNERGY SIII 2KVA


8. Printer Brother HL54440D, 257613
To print sample analysis results obtained by the instrument

9. KIT STACKING PHX50, 443910:

- Bracket – 2 pcs
- Washer – 12 pcs
- Screw – 12 pcs
- Leveling foot – 2 pcs
- Drawing
- Anchor – 2 pcs
- Nut hex – 2 pcs


to secure two Phoenix M50 instruments together in a stacked configuration

10. SEISMIC KIT PHX50, 443918:

- Bracket – 5 pcs
- Screw – 12 pcs
- Washer plain M5.3 X 10 X 1.0, SS – 12 pcs
- Bolt hex – 8 pcs
- Nut hex – 4 pcs
- Washer plain 656X 11.31 X.095 – 8 pcs
- Washer lock HSR. 625 – 8 pcs


to secure the Phoenix M50 instrument to a laboratory bench so it will not fall during an earth quake

11. ADAPTER PLATE FX40 TO PHX50 SPARE, 443873:

- Bracket – 6 pcs
- Screw – 12 pcs
- Washer – 8 pcs

to secure a BACTEC FX40 blood culture instrument to the top of a Phoenix M50

15. Peculiarities of exploitation

Each day it is necessary to perform routine maintenance procedures for BD Phoenix 100. The most suitable time for routine maintenance is morning, but these procedures may be performed in any convenient time. It is necessary to check the following elements:

1. Presence of paper in the printer. If it is few paper or it is lacking, substitute paper, following manufacturer's exploitation instructions.
2. Temperature readings on LCD display and temperature standard panel.

Please keep in mind that temperature panel readings may be seen by selection of one of the functions of device indicator control in technical maintenance menu. The temperature should lie within the range of 35 ± 1.5 °C.

The BD Phoenix M50 instrument requires little maintenance from the user to provide reliable performance. Daily activities include checking the instrument temperature and printer paper supply. On a weekly basis check the operation of the station status indicators (LEDs), the audible alarm, and the alert indicator. Routine preventive maintenance consists of a monthly check of the ambient air filter. All other procedures are on an as needed basis. Any maintenance or repair not described in this section should be performed by BD personnel only.

16. Safety measures

- a. It is necessary to ensure free access to the intake filter, located in the upper part of Phoenix device. Limitation of air flow may cause raise of the device temperature, which may affect the results of testing and cause equipment damage.
- b. Do not try to raise the device without using mechanical loading unit, as the device has big mass and dimensions. The attempt to raise the device without the appropriate facilities may cause serious traumas.
- c. Before start working with the device, all system users should attentively get acquainted with all control elements and indicators.
- d. Keep the prescribed safety measures, aimed against microbiological threats, when performing these procedures. All samples should be processed in accordance with CDC-NIH recommendations, SIN instructions or with the instructions of the institution on work with potentially infected human serum, blood or other human liquids. Before recycling, sterilize containers for samples and other infected materials using autoclaving.
- e. Besides gloves, it is recommended to use disposable laboratory overalls and protective or safety glasses when working with the device.
- f. The door of the device is blocked with software controlled electromechanic lock. Do not try to switch off door locking mechanism or to open the door, if "unlocked" sign is not seen. Rotating slewing mechanism may cause personal injury. If the slewing mechanism is not stopped when opening the door, immediately contact BD company for maintenance. Do not try to manually rotate the slewing mechanism; it may cause personal injury.
- g. When the system notifies you about the mistake or warning, it is necessary to immediately take the respective steps.
- h. Technical maintenance procedures or repair should be performed by qualified specialists. Failure to observe the present warning may entail accident or device irregularity.
- i. When installing the old filter, be sure that the filter is in the same position, as during its extraction (for example, the arrow on the filter points in the same direction).
- j. Before disinfection start, be sure that all parts of the body, which may come in contact with the engaged device surfaces, are completely covered.

k. In case of appearance of any subcode, not resulting in the list below, please put it down and contact BD company.

17. Special instructions

For diagnostic *in vitro* use.
All patients' samples and microbial cultures arte potentially infected and require observance of universal safety measures. Consult CDC Bio-safety in Microbiological and Biomedical Laboratories (Biological safety in microbiological and biomedical laboratories), 4th edition, 1999, as well as other recommended literature.
After seeding, one should be attentive with panels until they are put into the device.

18. Maintenance and current repair

Reliable operation of Phoenix system does not require large amounts of technical maintenance. Daily control of device temperature and presence of paper for printer is necessary. Each week it is necessary to check operation of station state indicators (LED), acoustic alarm and system warning indicator. Standard preventive technical maintenance includes monthly check of circulation system air filter. All other procedures should be performed as and when necessary.

For detailed information about maintenance rules see clause 6 ("Technical maintenance") of the Manual for BD Phoenix 100 and clause 5 (Maintenance) of the Manual for BD Phoenix M50. Technical maintenance or repair not described in the present sections of the Manual should be performed only by BD personnel.

19. Packing

Transport packing:
Box made of corrugated cardboard with inserts - transport protection devices (distance bars made of foam blocks).

20. Marking and symbols in technical documentation

Symbols	
	Manufacturer
	European representative
	Code number
	Series number
	CE mark
	Dangerous! Observe such warnings as obligatory condition for traumatism prevention. Read Reference Guide
	Electric contacts with dangerous voltage level!
	Temperature mode
	For <i>In Vitro</i> diagnostics
	Read the instruction

	On / Off
	Biological hazard
	Connection for earth wire

21. Transportation

The microbiology system should be transported and carried only in the original packing, taking into account rules of transportation in force for the given mode of transport.

22. Storage

Storage conditions

Temperature: 17.8 – >65.0 °C

Humidity: Relative humidity: 20-80% without condensing.

23. Operating life

By evaluation of the service history of the BD Phoenix 100 instrument, the useful life of the BD Phoenix M50 is 7 years. Through maintenance, repair and replacement, the BD Phoenix M50 could be used much longer than 7 years.

24. Safety requirements when working with the products

The system should be used only by specially prepared personnel in observance of safety rules for work with potentially infected material (CDC Bio-safety in Microbiological and Biomedical Laboratories Guide (Biological safety in microbiological and biomedical laboratories, 4th edition, 1999, or similar recommended literature, not conflicting with state and local normative acts).

Do not use the device in case of failures, appropriate warning technical messages, defects of power-carrying parts insulation, as well as network cable, and case ruptures.

25. Requirements to environment protection in course of products exploitation

The product during shelf life does not have any environmental impact.

26. Rules of recycling

The owner / operator of BD Phoenix automated microbiology system is liable for recycling of all wastes (sample materials), accumulated in course of measurements, in compliance with legislative and local norms, as well as with protocols of the medical institution.

Recycling is performed in accordance with the active legislative provisions after expiration of BD Phoenix automated microbiology system shelf life as electronic equipment wastes.

It is recommended to recycle the system in compliance with Public health regulations (SanPiN) 2.1.7.2790-10 "Sanitary and epidemiological requirements to medical wastes treatment".

27. Limited manufacturer's warranty

This warranty represents special legal rights. Besides that, you may have other rights, differing from region to region.

The primary purchaser of BD Phoenix 100 automated microbiologic system is guaranteed absence of defects in materials and manufacture within one year after installation. Exclusive responsibility of BD company under the above warranty also foresees repair and substitution of devices and their components (save such consumables as cartridges for printers, paper or filters), which in normal operation conditions have become defective within one year after supply.

BD company will supply you new or repaired components at its own discretion. All components, supplied for substitution, will correspond to new constructive parameters and will have guarantee until expiration of one-year term. Substituted components become the property of BD company. It is necessary that the equipment subject to the present agreement should be installed in compliance with the recommendations and instructions of the Reference Guide of BD Phoenix system.

Faults of BD Phoenix system, caused by connection or disconnection of cables, coupling the device with systems, not agreed or not supplied by BD company, or by failure to observe due care when working with the device and system maintenance, entail termination of the warranty and discharge the manufacturer from the obligations under the present document.

The present warranty substitutes all other expressed or implied warranties, including marketability warranty or suitability for certain use. BD company under no circumstances is liable for indirect, incidental, actual or consequential losses, irrespective of presence or absence of notification by BD company about such losses.

The BD Phoenix M50 Automated Microbiology System is warranted to the original purchaser to be free from defects in materials and workmanship for a period of one year following installation. BD's sole responsibility under this warranty shall be to repair or replace any instrument or its components (except for expendable supplies such as printer cartridges, paper, or filters) which under normal operating conditions, prove to be defective within one year of delivery.

BD will furnish new or re-manufactured components upon its option. All replacements shall meet new part specifications and shall be warranted as above for the remainder of the one year period. Replaced components become the property of BD.

It is understood that the equipment covered by this Agreement has been installed in accordance with the recommendations and instructions in the BD Phoenix M50 System User's Manual.

Any damage to a BD Phoenix M50 instrument resulting from the insertion or removal of cables that connect this instrument to systems other than those approved or supplied by BD, or the failure of the owner to maintain reasonable care and precautions in the operation and maintenance of the system, will void this warranty and terminate the obligations of the manufacturer as stated herein.

This warranty is in lieu of all other warranties, whether express or implied, including but not limited to, warranties of merchantability, or fitness for a particular use. In no event will BD be liable for indirect, incidental, special or consequential damages regardless of whether BD has been advised of such.

28. Information about development and manufacture of the medical product

BD Phoenix Automated Microbiology System, variations: BD Phoenix 100, BD Phoenix M50, with Accessories, manufacturer:

Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152, USA, should correspond to the internal technical regulation of Becton Dickinson and Company, USA (being confidential commercial information as containing technological know-how) on development and manufacture of BD Phoenix Automated Microbiology System.

The right for development of the medical product belongs to Becton Dickinson and Company, USA.

29. Denomination and legal address of the authorized representative

"Becton Dickinson Vostok" Limited Liability Company.

("Becton Dickinson Vostok" LLC)

Legal address: Block C, bldg. 1, 12, ul. Dvintsev, Moscow, Russian Federation, 127018

Actual address: Block C, bldg. 1, 12, ul. Dvintsev, Moscow, Russian Federation, 127018

Tel/fax: +7 495 775 85 82

«APPROVE»

Becton Dickinson and Company,
7 Loveton Circle, Sparks, Maryland 21152, USA

Date: *November 28, 2016*

Carol C. Tedford Name
(Signature)

Seal

Jennifer A. Creamer
11/28/2016

Jennifer A. Creamer
Notary Public, Maryland
Baltimore County
My Commission Expires:
October 19, 2019

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Becton Dickinson and Company,
7 Loveton Circle, Sparks, Maryland 21152, USA

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Becton Dickinson and Company, 7 Loveton Circle, Sparks, Maryland 21152,
USA

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11/30/07 09:11
VX.XXY / VA.AAB

Instrument # 1
Serial # 12345

Instrument Temperature Pass
Carousel Rotational Test Pass
Power Supply Check Pass

Normalizer Panels

<u>Tier</u>	<u>Sequence #</u>	<u>Status</u>	<u>Expiration Status</u>
A	429990000001	Pass	10/21/08
B	429990000002	Pass	10/21/08
C	429990000003	Pass	10/21/08
D	429990000004	Pass	10/21/08

Daily

	<u>Pass</u>	<u>Fail</u>	<u>Tech ID</u>
Record Instrument Temperature From Main Screen Range: (33.5 - 36.5)	_____	_____	_____
Record Standard Panel Temperature Range: (33.5 - 36.5)	_____	_____	_____
Check Printer Paper Supply	_____		

Weekly

Instrument Internal Green LEDs	_____	_____	_____
Instrument Internal Red LEDs	_____	_____	_____
Instrument Internal Amber LEDs	_____	_____	_____
Instrument System Alert Indicator	_____	_____	_____
Instrument Audible Alarm	_____	_____	_____

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BDXpert Rule Set Database Report			
		03/17/05 09:11	
Rule Set Based On	CLSI		
<u>Rule #</u>			
101	Only results for ampicillin, a fluoroquinolone, and trimethoprim-sulfamethoxazole should be reported for fecal isolates of <i>Salmonella</i> spp. or <i>Shigella</i> spp. No other drugs should be reported.	Enabled	Automatic
103	Results for chloramphenicol and a third generation cephalosporin should be reported for <i>Salmonella</i> spp. from extraintestinal sites.	Enabled	Automatic
104	<i>Salmonella</i> spp. and <i>Shigella</i> spp. should not be reported as susceptible or intermediate to first and second generation cephalosporins even if they appear to be susceptible in vitro because they are not effective clinically.	Enabled	Automatic
105	Ampicillin is the class representative for amoxicillin.	Enabled	Automatic
106	-		
.			
.			
.			
End of Report			

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- 1 Bronfenbrenner, J., and M.J. Schlesinger. 1918. A rapid method for the identification of bacteria fermenting carbohydrates. Am. J. Public. Health. 8:922-923.
- 2 Arnold, W.M., Jr., and R.H. Weaver. 1948. Quick microtechniques for the identification of cultures. I. Indole production. J. Lab. Clin. Med. 33:195-195.
- 3 Bachmann, B., and R.H. Weaver. 1951. Rapid microtechniques for identification of cultures. V. Reduction of nitrates. Am. J. Clin. Pathol. 21:195-196.
- 4 Hannan, J., and R.H. Weaver. 1949. Quick microtechniques for the identification of cultures II. Germentations. J. Lab. Clin. Med. 33:1338-1341.
- 5 Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
- 6 Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
- 7 Snyder, M.L. 1954. Paper discs containing entire culture medium for the differentiation of bacteria. Pathol. Bacteriol. 67:217-226.
- 8 Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Publ. Hlth. Trop. Med. 25:96-100.
- 9 Weaver, R.H. 1954. Quicker bacteriological results. Am. J. Med. Technol. 20:14-26.
- 10 Kampf, P., O. Rauhoff, and W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877-2879.
- 11 Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
- 12 Rammelkamp, C.H., and T. Maxon. 1942. Resistance of *Staphylococcus aureus* to the action of penicillin. Proc. Soc. Biol. And Med. 51:386-389.
- 13 Marymont, J.H. and R.M. Wentz. 1966. Serial dilution antibiotic sensitivity testing with the microtiter system. Am. J. Clin. Pathol. 45:548-551.
- 14 Gavan, R.L., and M.A. Town. 1970. A microdilution method for antibiotic susceptibility testing: an evaluation. Am. J. Clin. Pathol. 53:880-885.
- 15 Lancaster, M.V. and D. Rebecca. 1996. Antibiotic and cytotoxic drug susceptibility assays using resazurin and poisoning agents. US Patent #5,501,959.

- 16 Данные доступны в компании BD Diagnostics
- 17 Murray, Patrick R., et al. ed., *Manual of Clinical Microbiology*, 7th Edition, ASM Press, Washington, D.C., 1999.
- 18 DM Livermore. 1995. Beta-Lactamases in laboratory and clinical resistance. *Clin. Microbiol. Rev.* 8: 557-584.
- 19 K Bush, GA Jacoby, and AA Medeiros. 1995. A functional classification scheme for beta-lactamases and its correlation with molecular structure. *Antimicrob. Agents Chemother.* 39: 1211-1233.
- 20 GA Jacoby and P Han. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. *J. Clin. Microbiol.* 34: 908-911.
- 21 Kenneth S. Thomson, Christine C. Sanders, and Ellen Smith Moland. 1999. Use of Microdilution Panels with and without β -Lactamase Inhibitors as a Phenotypic Test for β -Lactamase Production among *Escherichia coli*, *Klebsiella* spp., *Enterobacter* spp., *Citrobacter freundii*, and *Serratia marcescens*. *Antimicrob. Agents Chemother.* 43: 1393-1400.
- 22 CLSI. M100-S18 *Performance Standards for Antimicrobial Susceptibility Testing; Eighteenth Informational Supplement*. January, 2008.
- 23 CLSI. M07-A8 *Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard-Eighth Edition*. January, 2009.
- 24 Livermore, D.M., T.G. Winstanley, and K.P. Shannon. 2001. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. *J. Antimicrob. Chemother.* 48, Suppl. S1, 87-102.
- 25 Cauwelier, et al. "Evaluation of a disk diffusion method with cefoxitin (30 mcg) for detection of methicillin-resistant *Staphylococcus aureus*." *Eur J. Clin. Microbiol Infect Dis.* May 2004,23(5):389-392.
26. Barnett, J.A., Payne, R.W. and Yarrow, D., *Yeasts: Characteristics and Identification, 2nd Edition*, Cambridge University Press, Cambridge, U.K. 1990.
27. Larone, D. H., *Medically Important Fungi, A Guide to Identification, 4th Edition*, ASM Press, Washington, D.C. 2002
28. de Hoog, G.S., Guarro, J., Gené, J. and Figueras, M.J., *Atlas of Clinical Fungi, 2nd Edition*, Centraalbureau voor Schimmelcultures, Utrecht, the Netherlands/Universita Rovira I Virgili, Reus, Spain, 2000.

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- 4 Hannan, J., and R.H. Weaver. 1949. Quick microtechniques for the identification of cultures II. Germentations. J. Lab. Clin. Med. 33:1338-1341.
- 5 Hartman, P.A. 1968. Miniaturized microbiological methods. Academic Press, New York.
- 6 Sanders, A.C., J.E. Faber, and T.M. Cook. 1957. A rapid method for the characterization of enteric pathogen using paper discs. Appl. Microbiol. 5:36-40.
- 7 Snyder, M.L. 1954. Paper discs containing entire culture medium for the differentiation of bacteria. Pathol. Bacteriol. 67:217-226.
- 8 Soto, O.B. 1949. Fermentation reactions with dried paper discs containing carbohydrate and indicator. Puerto Rican J. Publ. Hlth. Trop. Med. 25:96-100.
- 9 Weaver, R.H. 1954. Quicker bacteriological results. Am. J. Med. Technol. 20:14-26.
- 10 Kampfer, P., O. Rauhoff, And W. Dott. 1991. Glycosidase profiles of members of the family *Enterobacteriaceae*. J. Clin. Microbiol. 29:2877-2879.
- 11 Manafi, M., W. Kneifel, and S. Bascomb. 1991. Fluorogenic and chromogenic substrates used in bacterial diagnostics. Microbiol. Rev. 55:335-348.
- 12 Rammelkamp, C.H., and T. Maxon. 1942. Resistance of *Staphylococcus aureus* to the action of penicillin. Proc. Soc. Biol. And Med. 51:386-389.
- 13 Marymont, J.H. and R.M. Wentz. 1966. Serial dilution antibiotic sensitivity testing with the microtiterator system. Am. J. Clin. Pathol. 45:548-551.
- 14 Gavan, R.L., and M.A. Town. 1970. A microdilution method for antibiotic susceptibility testing: an evaluation. Am. J. Clin. Pathol. 53:880-885.
- 15 Lancaster, M.V. and D. Rebecca. 1996. Antibiotic and cytotoxic drug susceptibility assays using resazurin and poisoning agents. US Patent #5,501,959.
- 16 Данные компании Becton, Dickinson and Company.
- 17 Jorgensen, James H., et al. ed., Manual of Clinical Microbiology, 11th Edition, ASM Press, Washington, D.C., 2015.
- 18 DM Livermore. 1995. Beta-Lactamases in laboratory and clinical resistance. Clin. Microbiol. Rev. 8: 557-584.
- 19 K Bush, GA Jacoby, and AA Medeiros. 1995. A functional classification scheme for beta-lactamases and its correlation with molecular structure. Antimicrob. Agents Chemother. 39: 1211-1233.
- 20 GA Jacoby and P Han. 1996. Detection of extended-spectrum beta-lactamases in clinical isolates of *Klebsiella pneumoniae* and *Escherichia coli*. J. Clin. Microbiol. 34: 908-911.
- 21 Kenneth S. Thomson, Christine C. Sanders, and Ellen Smith Moland. 1999. Use of Microdilution Panels with and without Beta-Lactamase Inhibitors as a Phenotypic Test for Beta- Lactamase Production among *Escherichia coli*, *Klebsiella* spp., *Enterobacter* spp., *Citrobacter freundii*, and *Serratia marcescens*. Antimicrob. Agents Chemother. 43: 1393-1400.
- 22 CLSI. M100-S26 Performance Standards for Antimicrobial Susceptibility Testing; Twenty-sixth Informational Supplement. January 2016.
- 23 CLSI. M07-A10 Methods for Dilution Antimicrobial Susceptibility Tests for Bacteria That Grow Aerobically; Approved Standard – Tenth Edition. January 2015.
- 24 Livermore, D.M., T.G. Winstanley, and K.P. Shannon. 2001. Interpretative reading: recognizing the unusual and inferring resistance mechanisms from resistance phenotypes. J. Antimicrob. Chemother. 48, Suppl. S1, 87-102.
- 25 Cauwelier, et al. "Evaluation of a disk diffusion method with cefoxitin (30 mcg) for detection of methicillinresistant *Staphylococcus aureus*." Eur J. Clin. Microbiol Infect Dis. May 2004,23(5):389-392.
- 26 Barnett, J.A., Payne, R.W. and Yarrow, D., Yeasts: Characteristics and Identification, 3rd Edition, Cambridge University Press, Cambridge, U.K. 2000.
- 27 Larone, D. H., Medically Important Fungi, A Guide to Identification, 5th Edition, ASM Press, Washington, D.C. 2011.
- 28 de Hoog, G.S., Guarro, J., Gené, J. and Figueras, M.J., Atlas of Clinical Fungi, 4th Edition, Centraalbureau voor Schimmelcultures, Utrecht, the Netherlands/Universita Rovira I Virgili, Reus, Spain, 2014.

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