MILESTONE
H E L P I N G
P A T I E N T S

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

High Vacuum Biospecimens Transfer System TissueSAFE plus, in versions, with accessories

Operator Manual – MM123

APPROVED

Milestone S.r.l., Italy

(CORTESI Diego, General Director)

MILESTONE S.r.l.
Via Fatebenefratelli, 1/5
24010 SORISOLE (BG)
P. IVA 01879330163

This Operator Manual is related to High Vacuum Biospecimens Transfer System TissueSAFE plus, in versions, with accessories (hereafter referred to as the unit, the device).

We would like to thank you for having chosen our system, and welcome you to the growing club of Milestone equipment users. We're certain that you'll be completely satisfied with this new instrument's performance at your laboratory.

Please read this operator manual carefully and always keep it on hand for quick consultation.

Please don't hesitate to contact our representative in your country for any clarifications or support you may require:

Authorized Representative in the Russian Federation: BioLine, LLC

Legal address: Russia, 197022, Saint Petersburg, inner-city Municipality Aptekarsky Ostrov,
23 Professora Popova Str., Let. E, premises 3H/5H, room 208

Postal address: 197022, Saint Petersburg, 23 Professora Popova Str., Let. E

tel.: +7(812) 320-49-49

e-mail: main@bioline.ru

Or contact:

Manufacturer: Milestone S.r.l.

address: Via Fatebenefratelli, 1/5, IT-24010 Sorisole (BG), Italy

Web: www.milestonesrl.com

tel.: +39 035 571067

fax: +39 035 575498

Manufacturing site: Milestone S.r.l.

Corso Europe Unita, 6

IT-24030 Valbrembo (BG), Italy

Read this manual carefully before using the instrument.

Summary

1. INTRODUCTION	5
1.1. Symbols used in this manual	5
1.2. Safety requirements and measures.....	5
1.3. General description.....	8
1.4. Structure and functioning	10
1.5. Technical specification	13
1.6. Touchscreen terminal	16
1.7. Warnings	16
1.8. Information on labels and packaging of the unit	18
1.9. Barcode reader.....	18
2. INSTALLING THE INSTRUMENT	19
2.1. Space requirements.....	19
2.2. Unpacking and inspection	20
2.3. How to insert the shelves into the holder for bags.....	21
2.4. How to install the barcode reader	22
2.5. Fume extraction system	23
2.5.1. Power supply	23
2.5.2. Rear panel connections.....	24
2.6. Final positioning.....	24
2.6.1. Activation/Deactivation.....	25
3. MILESTONE VACUUM BAGS	26
4. USING THE TISSUESAFE PLUS UNIT	27
4.1. How to login and logout from the system.....	27
4.1.1. Creating a new user	28
4.2. Setting the date and time	29
4.3. Owner settings.....	30
4.4. Executing a default program.....	34
4.4.1. Using one bag more than once for the same patient (single use).....	41
4.4.2. Power failure during the program	41
4.5. How to set favourite programs	42
4.6. Recovering the results of logged processes	43
4.6.1. How to delete saved runs.....	45
4.6.2. Exporting data from the terminal.....	45
4.6.3. Installing and using LogViewer.....	47
4.7. How to modify a program	49
4.8. How to set favourite programs	50
4.9. Temperature Monitoring System.....	51
4.10. Sample transportation.....	56
4.11. How to use the bag sealer.....	57
5. COMPONENTS AND ACCESSORIES	57

6. ADDITIONAL DEVICES THAT ARE NOT PART OF THE INSTRUMENT, BUT INTENDED FOR CO-USE	63
6.1. Milestone Vacuum Bags.....	63
6.1.1. Storing the vacuum bags.....	64
6.1.2. Stand up vacuum bags	64
6.1.3. Vacuum bags (not stand up)	64
6.2. Temperature indicators.....	65
7. MAINTENANCE.....	66
7.1. Operating room personnel	66
7.1.1. Daily	66
7.1.2. Weekly.....	66
7.2. Service staff.....	67
8. WASTE DISPOSAL	68
8.1. Disposal of vacuum bags.....	68
8.2. Waste of the unit, its components and accessories	68
9. DECONTAMINATION AND CLEANING	69
9.1. Safety precautions and regulations	69
9.2. Procedure for decontaminating and cleaning the TissueSAFE plus unit.....	70
9.2.1. Detergents.....	70
9.2.2. Procedure for decontaminating and cleaning the unit.....	70
10. TRANSPORTATION AND STORAGE	72
11. MANUFACTURER WARRANTY	73
A. APPENDIX A – PROPER OPERATING METHODS.....	74
A.2. HOW TO INSERT A SAMPLE INTO A VACUUM BAG.....	74
A.3. HOW TO PROPERLY SEAL A VACUUM BAG	74
B. APPENDIX B – VACUUM PROGRAMS GUIDELINES	75
C. APPENDIX C – TRANSPORTATION GUIDELINES	76
D. APPENDIX D – PRINTER AND BARCODE READER	77
D.1. How to use the printer	77
D.1.1. Storing the rolls of paper	78
D.1.2. Changing the roll of paper	78
D.1.3. Printer troubleshooting	79

1. INTRODUCTION

1.1. Symbols used in this manual

	An instruction accompanied by this symbol provides a caution against electric shock. Be sure to follow the instruction to prevent accidents resulting from electric shock.
	Biohazard: be careful when you execute the procedure marked with this symbol: danger of biohazard contamination.
	Surfaces identified with this symbol are at danger temperature (more than 70°C).
	Laser light – Do not stand in front of the laser beam – CLASS II Laser Device.
	CE Logo: this instrument is compliant with the European directives.
	WEEE European directive symbol (2012/19/EU): Electric/electronic device do not throw out in the environment.
	IVD medical device according to 98/79/EC directive.
	Indicates the manufacturer's batch number.
	Indicates the manufacturer's serial number.
	Indicates the manufacturer's catalogue code.
	Fragile
	Manufacturer
	Manufacture date
	Storage upper limit of temperature.
	Storage humidity limits.
	Keep away from sunlight.
	See the Operator manual
	Best before
	Warning! Read and follow these rules carefully.
	Information with this symbol means that this function can only be performed by authorized and qualified TECHNICAL personnel

1.2 Safety requirements and measures

 Warning
➤ The safety and caution notes in this chapter must be observed at all times. Please read it, even if you are aware of the operation and use of other Milestone devices. ➤ It is forbidden to remove or modify the protective devices on the unit, its components and accessories. ➤ Only service personnel qualified by Milestone S.r.l. may repair the instrument and access the instrument's internal components. ➤ Residual risks: ➤ the unit is designed in accordance with the state of the art and recognized safety regulations. Improper use and handling may endanger the health and life of the user or third parties, as well as this can lead to damage to the device and property. ➤ The unit may only be operated for its intended use and in perfect working order for safety reasons. ➤ In the event of failures that could impair the level of safety, immediately take the unit out of service and contact the Milestone Customer Support Service. ➤ Only original spare parts and authorized Milestone accessories may be used. ➤ The requirements to electromagnetic compatibility, harmful emissions and immunity to interference, as well as the requirements in accordance with IEC 61326-2-6, must be met.

This unit is manufactured and tested in accordance with the safety regulations for electrical equipment for measurement, control, and laboratory use. To maintain this condition and ensure safe operation, the user must follow all instructions and warnings.

 Warning
➤ The unit is intended only for connection to a network with protective grounding. Milestone S.r.l. is not responsible for faults caused by connection to an unsecured network. ➤ Only technicians trained and authorized by Milestone are allowed to connect USB input devices. This also applies to the network connection, which can only be used in conjunction with remote diagnostics (service diagnostics).

 Warning
CLASS II Laser Device Laser light! Do not look at the beam! This can lead to retinal injury.

 Warning
 LASER LIGHT - DO NOT LOOK AT THE BEAM $P < 1 \text{ mWt}$, $\lambda = \text{up to } 670 \text{ nm}$ Pulse width = 12 mcs CLASS II Laser Device

Cautions

The protective devices fitted by the manufacturer to this unit are only the basis for the protection against accidents. The main responsibility for a fault-free working process shall be borne primarily by the owner of the institution where the device is operated, as well as the persons designated by it who operate, maintain or repair the unit.

To ensure perfect operation of the unit, the following instructions and warnings must be observed.

Please be aware that direct or indirect contact with the TissueSAFE plus can result in electrostatic discharge.

 Warning
During normal use, some surfaces of the device become hot. These surfaces are marked with warning symbols. Touching these surfaces without suitable safety measures can cause burns.

Warnings - Transportation and installation

 Warning
➤ The unit may be transported only in an upright position. ➤ Weight of the unit is 135 kg. ➤ Milestone S.r.l. recommends to involve a specialized transport company to transport, install or possible handling of the unit. ➤ Keep the packaging of the unit. ➤ Do not expose the device to direct sunlight. ➤ The unit may only be connected to a grounded outlet. Do not use extension cords without a protective conductor. ➤ Condensation water may form in the instrument if there is an extreme difference in temperature between the storage location and the installation site and if air humidity is high at the same time. If this is the case, wait at least two hours before switching on the instrument. ➤ The installation of the unit on site must only be carried out in the presence and under the supervision of a person trained by Milestone S.r.l. This also applies to the potential transport to a new area of use. ➤ It is recommended to involve a person trained by Milestone S.r.l. for re-commissioning of the unit.

Danger warnings - operation of the unit

 Warning
➤ The instrument may be operated by trained laboratory personnel only. It must only be operated for the purpose of its designated use and according to the instructions contained in this Operator manual. Wear antistatic protective clothing (e.g. out of natural fabrics) when operating the device. ➤ Wear suitable protective clothing (gown, goggles and gloves) when operating the unit to protect against contact with potentially infectious microbiological materials. ➤ In the event of an abnormal condition, turn off the power switch and disconnect the unit from the power supply. ➤ In the event of a serious fault of the unit, follow the warnings and error messages shown on the display. ➤ Make sure the unit is equipped with a carbon filter. In addition, we strongly recommend to connect the unit to a fume extraction system.

 Note	Milestone S.r.l. recommends to provide for an airflow of at least 12m ³ /h (cubic meters per hour) in the laboratory.
---	--

 Warning
➤ During operation the cover of the device must be closed. Milestone S.r.l. is not responsible for the deterioration in the quality of the unit, caused by opening the cover during the process of packaging and sealing biological samples in bags under vacuum. ➤ CAUTION when closing the cover: Risk of crushing! Do not put your hands in the area of cover movement. ➤ During operation and when cleaning the unit, liquid must not get into the cover or into the gaps.

Warnings - Maintenance and Cleaning

 Warning
➤ Switch off the unit and disconnect the plug from the power outlet before performing maintenance. ➤ Wear suitable protective clothing (lab coat and gloves) when cleaning the unit. ➤ Do not remove the side panels of the instrument during cleaning procedures. ➤ Disconnect the device from its electrical power supply before performing any operations. Do not use direct or high-pressure water to clean the instrument. ➤ DO NOT use alcohol-based detergents to clean the cover. ➤ The outer surfaces of the unit housing, components and accessories are suitable for disinfection with a 3% hydrogen peroxide solution with 0.5% universal detergent.

Safety devices fitted in the unit

The device can be moved only in the horizontal plane (X and Y axes).

 Warning
➤ The device can be moved only in the horizontal plane (X and Y axes). ➤ The unit comes with replaceable fuses, which must only be replaced by authorized technical personnel. If the system is not functioning properly, please contact the local representative of the manufacturer. All of the instrument's parts and accessories must be supplied exclusively by the manufacturer.

1.3. General description

Name of the medical device: High Vacuum Biospecimens Transfer System TissueSAFE plus, in versions, with accessories

Versions:

I. High Vacuum Biospecimens Transfer System TissueSAFE plus, basic instrument, including:

1. High Vacuum Biospecimens Transfer System TissueSAFE plus, basic instrument - 1 pc.;
2. Discharge tube - 1 pc.;
3. Metal strap - 2 pcs;
4. Main power cable - 1 pc.;
5. Bar code reader mounted on its support - 1 pc.;
6. Plastic spacers positioned in the cavity - 2 pcs;
7. Shelves to be inserted into the holder for bags - 3 pcs;
8. USB flash drive containing "LogVIEWER" software - 1 pc.
9. Rolls of printer paper - 10 pcs;
10. Stylus for working with the touch screen terminal - 1 pc.;
11. Operator Manual in electronic format - 1 pc.;
12. Bag sealer with handle - 12 pcs;
13. Operator Manual in paper format - 1 pc.;
14. Kit Refrigerator box in PS, 32 liters complete with N.4 cool packs - 1 pc.;

Accessories:

1. HEPA filter - 1 pc.;
2. Carbon filter - 1 pc.

II. High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system, including:

1. High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system - 1 pc.;
2. Discharge tube - 1 pc.;
3. Metal strap - 2 pcs;
4. Main power cable - 1 pc.;
5. Bar code reader mounted on its support - 1 pc.;
6. Plastic spacers positioned in the cavity - 2 pcs;
7. Shelves to be inserted into the holder for bags - 3 pcs;
8. USB flash drive containing "LogVIEWER" software - 1 pc.
9. Rolls of printer paper (if necessary) - 10 pcs;
10. Stylus for working with the touch screen terminal - 1 pc.;
11. Operator Manual in electronic format - 1 pc.;
12. Bag sealer with handle - 12 pcs;
13. Operator Manual in paper format - 1 pc.;
14. Kit Refrigerator box in PS, 32 liters complete with N.4 cool packs - 1 pc.;
15. Temperature monitoring system - 1 pc. (if necessary)

Accessories:

1. HEPA filter - 1 pc.;
2. Carbon filter - 1 pc.;

Version	CODE
High Vacuum Biospecimens Transfer System TissueSAFE plus, basic instrument	68900
High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system	68902

The instrument is only intended for laboratory use (indoors).
The TissueSAFE plus must only be used as described in this manual.

Any use of the instrument other than that which is permitted shall be considered improper and may void the product's warranty.

Intended use: TissueSAFE plus is an auxiliary equipment item for in vitro diagnostics, a specially designed vacuum system designated for packaging and sealing of biological samples and human organs in vacuum bags for closing of biological specimens on TISSUE SAFE and SealSAFE to maintain their stable state at +2...+4 °C for up to 72 hours for subsequent transportation and / or storage for further histological studies in histopathological laboratories of medical institutions.

Indications for use: For packaging and sealing of biological samples in bags under vacuum to maintain their stable state at +2...+4 °C for up to 72 hours for subsequent transportation and / or storage for further histological studies in histopathological laboratories of medical institutions.

Contraindications: There are no contraindications, if the user observes all the instructions and recommendations specified in the Operator manual.

Warnings: The unit is intended only for in vitro diagnostics. It is necessary to observe all the precautions normally required when handling and disposing of biological samples. Use disposable gloves when handling biological samples as they can be infected. Close vacuum bags containing biological samples carefully. Do not eat, drink or smoke while using the unit. The working surfaces of the unit must be cleaned regularly and thoroughly after each use. When installing the unit, it is necessary to take into account the electromagnetic compatibility (EMC) conditions and install it in strict accordance with the EMC information given in the Operator manual.

Potential side effects: not applicable

Information on potential consumers and demographical and population aspects of application: The unit is intended for use only by qualified specialists (pathologist, nurses, medical laboratory technicians). Use of the unit is possible both in large laboratories, hospitals, medical institutions with a large flow of tests, and in organizations with an average flow of tests. **Before starting operation, all users should be familiar with the Operator manual and trained to operate the unit.**

Type of the sample studied: biological samples of patient
Types of biological samples with which the unit is used: Type of samples - biological tissues and organic fluids

General information on biological samples used

List of biological sample types:

Kidney	Open colon	Stomach / esophagus
Breast	Pancreas	Ovary
Uterus	Bladder	Tonsil
Liver	Heart	Lung
Muscle	Placenta	Brain
Thyroid	Prostate	Lymph node
Thymus	Intestine/intestines	Closed colon
Skin	Testicle	Gallbladder
Tongue	Appendix	Eye

Spleen	Bone (not sharp)	Cyst <10cm in diameter
Scraping / TURP	Cyst >10cm in diameter	Fetus

Before vacuum packaging with TissueSAFE plus, the said biological samples, depending on their size, are placed (1 pc. per bag) in a disposable Milestone vacuum bag of appropriate size. Detailed information on how to use vacuum bags is provided in cl. 6. *"ADDITIONAL DEVICES THAT ARE NOT PART OF THE INSTRUMENT, BUT PROVIDED FOR CO-USE " of the Operation Manual - MM123.*

Depending on the type of histological study, the sample/organ can be placed in a vacuum bag either entirely or partially. Samples/organs containing biological tissues, mucous membranes, liquids, must be freely placed, taking into account their sizes, in a vacuum bag.

The said biological samples can be collected both from a living person during the surgery, and from the body of the deceased person on autopsy.

Place and conditions for the selection of biological samples, under which there is a need for their storage / transportation:

1. Operating Room. When a whole organ, its part, a cyst, etc. is extirpated during the surgery, the extirpated samples are packed in vacuum bags for subsequent transportation to a histological laboratory for further histological studies. If the surgery is performed on weekends or at night, it becomes necessary to store the sample until the opening of the histological laboratory.

2. Pathology Department. During the autopsy, samples of biological tissues of organs from the body are collected for subsequent histological studies in a histological laboratory. Samples are packed in vacuum bags for transportation and/or storage and subsequent transportation.

Stability of the biological samples after vacuum packaging with TissueSAFE plus: storage up to 72 hours at +2°C ... +4 °C.

1.4 Structure and functioning

Purpose of use: TissueSAFE plus is an auxiliary equipment item for in vitro diagnostics, a specially designed vacuum system designated for packaging and sealing of biological samples and human organs in vacuum bags for closing of biological specimens on TISSUE SAFE and SealSAFE to maintain their stable state at +2°C...+4 °C for up to 72 hours for subsequent transportation and / or storage for further histological studies in histopathological laboratories of medical institutions.

Principles of operation


A qualified IVD specialist prepares a biological sample - puts it in a disposable vacuum bag of appropriate size. Patient's information is entered by the specialist manually or by scanning a barcode on a referral. For this action the unit is equipped with a barcode reader. The operator selects the relevant packaging protocol on the screen of the unit and puts a bag with biological sample in the vacuum chamber of TissueSAFE plus. After closing the cover of the unit, the program of packaging starts automatically. Then, the vacuum pump begins to create a vacuum inside the chamber, and the bag is automatically sealed with integrated welding bar. The unit is equipped with a printer and after sealing the bag, the printer prints a label with the appropriate identification information of the biological sample, which [label] the operator places in a special pocket on the bag. Now the biological sample is ready for subsequent transportation and / or storage.

Product configuration: The instrument includes a stainless steel chamber that allows a vacuum to be generated by means of an internal pump. The welding bar seals and closes a special bag for storing biological samples under vacuum. The instrument is controlled by an electronic board and a touchscreen terminal. The terminal identifies and records all processes and manages the correct sealing procedure.


The main elements of the instrument:

- vacuum chamber,
- barcode reader,
- welding bar,
- printer for printing labels,
- display for managing,
- shelves for storing a stock of vacuum bags,
- slot for temperature monitoring system*

*(only for High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system)


Front view

1. Touch screen terminal
2. Stylus for working with the touchscreen terminal
3. Printer
4. Barcode reader
5. Silicon bar
6. Cover Handle
7. Spacer
8. Welding bar
9. Vacuum chamber
10. Cover open-close switch

Rear view

1. Power connector
2. Place for marking
3. Gas supply hole - not used
4. ETH/LAN connection
5. Connection for fume extraction system (for connecting the instrument to the laboratory's fume extraction system (if needed)).

Design of the bag sealer device

There is a welding bar with pneumatic cylinders in the chamber for sealing vacuum bags. At the end of the evacuation cycle, the welding bar is pressed against the silicone bar by means of pneumatic cylinders for a predetermined time, set by the device program, necessary for sealing the bag.

By means of a special transformer (230 V to 18 V), the electric current is converted into heat, the welding bar is heated and the unit seals the bag. The power is supplied for only a few seconds, which is controlled by the electronic board, in addition, the hardware control of the timer provides a safety control that turns off the power supply to the transformer in the event of an electronic board failure.

Spacers that regulate the depth of the chamber can be placed in the working area of the chamber (for sealing the vacuum bags of small and medium sizes).


The bag must be positioned in the chamber so that the welding bar is between two white lines printed on the bag, as shown in the figure below.

Always position the uppermost line over the welding bar.

If the centreline of the sample in bag does not coincide with the height of the welding bar, remove one or both of the spacers from the TissueSAFE plus unit's chamber.

Also, the camera has a plexiglass cover, inside of which there is a silicone bar (marked with A), for tightly pressing the bag to the welding bar. For a tighter fit of the cover to the steel surface, a rubber gasket is mounted around the perimeter (marked with B).

1.5. Technical specification

Standards applied

Directives applied:

- Regulation 2017/746 of the European Parliament and of the Council of 5 April 2017 on in vitro diagnostic medical devices
- Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment
- Directive 2014/35/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits
- Directive 2014/30/EU of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to electromagnetic compatibility

International standards applied

- EN ISO 14971 Medical devices -- Application of risk management to medical devices
- EN ISO 18113-1 – In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements;
- EN ISO 18113-2 – In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 2. In vitro diagnostic reagents for professional use;
- EN ISO 18113-3 In vitro diagnostic medical devices -- Information supplied by the manufacturer (labelling) – Part 3: In vitro diagnostic equipment for professional use
- EN ISO 15223-1 Symbols for use in the labeling of medical devices
- EN 61010-2-101 Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 2-101. Particular requirements for in vitro diagnostic (IVD) medical equipment
- EN 61326-2-6 – Electrical equipment for measurement, control and laboratory use - EMC requirements - Part 2-6. Particular requirements - In vitro diagnostic (IVD) medical equipment;

8. EN 62304 – Medical devices. Software. Life cycle processes;
9. EN 62366 – Medical devices - Application of usability engineering to medical devices.

Low Voltage Directive

CEI EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use.

Electromagnetic compatibility

1. CEI EN 61326-1 – Electrical equipment for measurement, control, and laboratory use:
 - a. CISPR 11 (CEI EN 55011);
 - b. EN 61000-3-2 – Harmonic current emission;
 - c. EN 61000-3-2 – Flicker;
 - d. EN 61000-4-2 – Electrical Static Discharges;
 - e. EN 61000-4-3 – Radio-frequency electromagnetic field;
 - f. EN 61000-4-4 – Electrical fast transient/burst;
 - g. EN 61000-4-5 – Surge;
 - h. EN 61000-4-6 – Conducted disturbances induced by RF fields;
 - i. EN 61000-4-11 – Voltage dips, short interruptions and voltage variations.

Russian standards

- GOST R 50444-2020 Medical instruments, apparatus and equipment. General specifications.
- GOST R ISO 18113-1-2015 – In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling).
- Part 1. Terms, definitions and general requirements
- GOST R ISO 18113-3-2015 In vitro diagnostic medical devices -- Information supplied by the manufacturer (labelling) –
- Part 3: In vitro diagnostic equipment for professional use
- GOST IEC 61010-1-2014 Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1. General requirements
- GOST IEC 61010-2-010-2013 Safety requirements for electrical equipment for measurement, control, and laboratory use.
- Part 2-010. Particular requirements for laboratory equipment for the heating of materials
- GOST IEC 60825-1-2013 Safety of laser products. Part 1. Equipment classification, requirements and user's guide
- GOST R IEC 62304-2013 Medical devices. Software. Life cycle processes
- GOST R IEC 62366-1-2021 Medical devices. Part 1. - Application of usability engineering to medical devices
- GOST R ISO 15223-1-2020 Medical devices. Symbols applied in medical device labelling, marking and in supporting documentation. Part 1. Essential Requirements
- GOST IEC 61326-1-2014 Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1. General requirements

Main parameters and technical specifications

Technical specifications of TissueSAFE plus, its components and accessories are given in table 1.

Table 1

Technical specification	TS value / parameter
Rated voltage:	220 V ~ +10 %
Rated frequency:	50 Hz
Max power consumption:	Not more than 1800VA
Fuses:	Dimensions (diameter x length) 10(±1) x 38 (±1) mm; Rated current: 8 A; Rated voltage: 250V; Color: white. Code: gG8A250V. Weight - 10 (±1) g. Quantity used - 2 pcs
Overvoltage category:	II
Maximum allowable time for starting the operation mode	Calculated from the moment of switching on or starting the device, it should be no more than: 45 seconds.
Continuous operation:	It is not recommended to carry out more than 3 cycles during 15 minutes to prevent thermal protection from tripping.
Mode of operation:	CONTINUOUS OPERATION
Characteristics of the vacuum chamber:	Evacuation: - Stainless steel chamber, vacuum resistant - EVC (Electronic Vacuum Control) System with electronic sensor for the best gas/vacuum characteristics in any atmospheric/altitude conditions - Volume of the vacuum chamber: 0,023 m ³ - Vacuum pump: vacuum pressure up to 10 mbar - Number of operation programs: 4: Program P1: Speed of vacuum creating (deflation) - 25 s starting from 927 mbar; Max. level of vacuum created in the chamber <u>25 mbar</u> ; Program P2: Speed of vacuum creating (deflation) - 19 s starting from 927 mbar; Max. level of vacuum created in the chamber <u>50 mbar</u> ; Program P3: Speed of vacuum creating (deflation) - 16 s starting from 927 mbar; Max. level of vacuum created in the chamber <u>75 mbar</u> ; Program P4: Speed of vacuum creating (deflation) - 10 s starting from 927 mbar; Max. level of vacuum created in the chamber <u>325 mbar</u> ; Sealing: - General description of the bag sealer device design is given above, in cl. 1.4 of this Operator Manual - Sealing seam width: double sealing line 4 mm - Speed of sealing: 3,3 s - melting temperature of film (BOPA / PE), from which the bags are made: 90 (± 10) ° C, Note! The unit does not control the sealing temperature. - requirements to the seam quality: no gaps in sealing, underheating, overheating
Slot for temperature monitoring system	Version: 1) High Vacuum Biospecimens Transfer System TissueSAFE plus - none 2) High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR - yes (Slot type: USB port 2.0)
Ingress protection rating as per GOST 14254	IP 20
Characteristics of the integrated printer	Information printed on the label: - Program P1-P4 (Name of the program) -Case ID (if present) -Bag ID -User (User who executed the program) -Date -Time -Outcome: Completed with success/Stopped & Sealed/Aborted. Type: Thermal printer Width of record: 58 mm Paper feed speed (speed of printing): 50 mm/s Type of record: Continuous real-time recording
Noise level:	not more than 70 dBA
Air flow through the fan:	not more than 12m ³ /h (cubic meters per hour)
Operation, transportation and storage conditions:	Operation: Temperature: +15°C... +30 °C. Relative humidity: 10-80%, non-condensing, altitude up to 2000 m Transportation and Storage: Temperature: from -20°C... +60 °C. Relative humidity: up to 50% (at 40 °C), non-condensing.
Dimensions (length x width x height):	With cover opened: 84 (±5%) x 57 (±5%) x 135 (±5%) cm With cover closed: 84 (±5%) x 57 (±5%) x 105 (±5%) cm

Weight	135(±10%) kg		
Fume extraction system	12 m ³ /h (7 cfm)		
Electromagnetic compatibility	The instrument meets the requirements of GOST R IEC 61326-1, GOST R IEC 61326-2-6		
	TS Classification		
	Reference document	Class	Group
	GOST R 51318.11	A	1
	GOST 30804.3.2	A	
	GOST 30804.4.11	1	

1.6. Touchscreen terminal

The following photo shows the TissueSAFE plus unit's touchscreen terminal. The unit can be controlled by simply touching the icons on the screen.

Two USB ports can be found on the sides of the instrument, one on the right and one on the left.

When the USB ports are not in use, the protective plug must be inserted into the USB slot, as shown in the following figure.

Characteristics:

TFT-LCD 7" colour display, 800X480 resolution, Resistive touchscreen, Mass storage 4Gb
 CPU 1GHz, OS used: Windows™ Embedded Compact 7

Main characteristics of the integrated software

Multicomponent software, customizable according to user needs

SW name	TissueSAFE plus
SW version	17.12 and higher
SW classification	Class A
release date	24.01.2018

I/O-Interface:

2 USB 2.0 ports, 1 LAN Ethernet 10/100 Mbps RJ-45

1.7. Warnings

The warranty will be void if the instrument is used for purposes other than those specified by the manufacturer. The plug is to be considered a disconnection device. Always disconnect the plug from the power outlet before assembling the system, connecting its accessories (where present), and performing cleaning operations. The device must be connected to the main power outlets, which must be visible and accessible to the operator.

The device comes with replaceable fuses, watch must only be replaced by authorized technical personnel. If the system is not functioning properly, please contact the local representative of the manufacturer.

All of the instrument's parts and accessories must be supplied exclusively by the manufacturer.

DO NOT OPEN THE CASING. Any failure to comply with this instruction could compromise the device's safety. If the device is not functioning properly, please contact the supplier or the manufacturer.

Do not remove the instrument, the accessories, the panels, or the covers. Any repairs and/or access to the instrument's internal components must only be carried out by qualified and authorized technicians.

It is important to comply with the safety regulations and the standards of good laboratory practices. It is recommended to always use common sense and the best known practices while using the instrument.

Together with the personnel in charge of its use, maintenance, and repair, the organization/institution that owns the instrument itself is primarily responsible for its use under safe conditions.

In order to avoid damaging the instrument, only use the vacuum bags described in this manual, which are supplied by Milestone s.r.l. For more details please refer to chapter 6.1. For disposal instructions please refer to chapter 8.1.

In order to avoid overheating the system, it is recommended to avoid performing more than three sample sealing operations every 15 minutes.

The welding bar can be extremely hot after welding operations.

Do not obstruct air discharge.

After being taken from its place of storage, the TissueSAFE plus must be acclimatized to the environmental conditions (at least half an hour) before being turned on.

1.8. Information on labels and packaging of the unit


The instrument's label

The following label can be found on the rear panel:

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The unit has an integrated printer that prints labels with the following information:

- Program P1-P4 (Name of the program)
- Case ID (if present)
- Bag ID
- User (User who executed the program)
- Date
- Time
- Outcome: Completed with Success/Stopped & Sealed/
Aborted.

For the symbols, please refer to chapter 1.1.

Information on packaging

The unit, its components and accessories are packed in polyethylene film, placed in a wooden crate of the following dimensions:

Width: 77cm ± 5%	Height: 134cm ± 5%	Depth: 102cm ± 5%
------------------	--------------------	-------------------

Weight in the packaging: 175 kg ± 10%.

The unit, its components and accessories in the crate are secured with bumper inserts that prevent free movement of the contents. Operational documentation is placed in a polyethylene bag and packed together with the unit. In such packaging, the unit is available for transportation and storage.

1.9. Barcode reader


The unit's package comes complete with a barcode reader.

This scanner is designed to allow for the quick and intuitive reading of barcodes. Its ergonomic design allows for hands-free use or manual operations.

Please refer to chapter 2.4 for more information on how to install it on the unit. The TissueSAFE plus unit's barcode reader is compatible with 1D and 2D codes.

2. INSTALLING THE INSTRUMENT

The TissueSAFE plus unit is a precision instrument that must be unpacked and installed with the utmost care. Move the TissueSAFE plus unit to its final destination. Make sure that the flooring is level and is not made from a flammable material. The following images show the dimensions of the TissueSAFE plus unit and the space requirements:

2.1. Space requirements

Environmental specifications:

Supply source: Supply source: Voltage: 230 V ($\pm 10\%$), 50 Hz AC

Power supply line needs to have a ground connection (Class I device).

Power line protection: 230V version: Miniature circuit breaker curve C, 10A, residual current device: 30mA.

Warning:

- Make sure the outlet is properly grounded. Improper grounding can result in electric shock and / or damage to the unit.
- Connect the unit only to a power outlet that meets the above requirements and has a proper fuse.
- Before connecting the plug to the laboratory's outlet, check to make sure that the unit's power supply voltage corresponds to that of the laboratory's network.
- The use of an incorrect power supply will result in serious damage.

The space requirements for positioning the unit are the following:

Width: 97cm $\pm 5\%$

Height: 155cm $\pm 5\%$

Depth: 145cm $\pm 5\%$

Weight of the unit: 135 kg $\pm 10\%$

Fume extraction system

The discharge of the TissueSAFE plus unit's vacuum pump, which is located on the back of the instrument, can be connected to a fume extraction system.

The supplied plastic tube has the following characteristics:

- Length: 3 m $\pm 5\%$
- External diameter: 58 mm $\pm 5\%$
- Internal diameter: 50 mm $\pm 5\%$

The external fume extraction system must allow for an airflow of at least 12m³/h (cubic meters per hour).

The unit comes with a HEPA filter and an active carbon filter, which are not accessible to the user. The first filter is designed to filter the vapours coming from the biological samples, while the second filter is designed to filter odours. The filters are changed every 6 months during the Preventive Maintenance intervention carried out by authorized technical personnel.

2.2. Unpacking and inspection

When unpacking the instrument it is necessary to verify that all the parts correspond to those which are listed on the shipping documents. The unit must be unpacked with the utmost care in order to avoid scratching or damaging its external casing.

The package includes the following elements:

I. High Vacuum Biospecimens Transfer System TissueSAFE plus, basic instrument, including:

1. High Vacuum Biospecimens Transfer System TissueSAFE plus, basic instrument - 1 pc.;
2. Discharge tube - 1 pc.;
3. Metal strap - 2 pcs;
4. Main power cable - 1 pc.;
5. Bar code reader mounted on its support - 1 pc.;
6. Plastic spacers positioned in the cavity - 2 pcs;
7. Shelves to be inserted into the holder for bags - 3 pcs;
8. USB flash drive containing "LogVIEWER" software - 1 pc.
9. Rolls of printer paper - 10 pcs;
10. Stylus for working with the touch screen terminal - 1 pc.;
11. Operator Manual in electronic format - 1 pc.;
12. Bag sealer with handle - 12 pcs;
13. Operator Manual in paper format - 1 pc.;
14. Kit Refrigerator box in PS, 32 liters complete with N.4 cool packs - 1 pc.;

Accessories:

1. HEPA filter - 1 pc.;
2. Carbon filter - 1 pc.

II. High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system, including:

1. High Vacuum Biospecimens Transfer System TissueSAFE plus WDLR, with integrated slot for temperature monitoring system - 1 pc.;
2. Discharge tube - 1 pc.;
3. Metal strap - 2 pcs;
4. Main power cable - 1 pc.;
5. Bar code reader mounted on its support - 1 pc.;
6. Plastic spacers positioned in the cavity - 2 pcs;
7. Shelves to be inserted into the holder for bags - 3 pcs;
8. USB flash drive containing "LogVIEWER" software - 1 pc.
9. Rolls of printer paper (if necessary) - 10 pcs;
10. Stylus for working with the touch screen terminal - 1 pc.;
11. Operator Manual in electronic format - 1 pc.;
12. Bag sealer with handle - 12 pcs;
13. Operator Manual in paper format - 1 pc.;
14. Kit Refrigerator box in PS, 32 liters complete with N.4 cool packs - 1 pc.;
15. Temperature monitoring system - 1 pc. (if necessary)

Accessories:

1. HEPA filter - 1 pc.;
2. Carbon filter - 1 pc.;

2.3. How to insert the shelves into the holder for bags

At the installation time it is required to insert the shelves into the holder for bags of the TissueSAFE plus.

In the figure below the empty holder is shown.

Insert the first shelf to be inserted into the holder for bags as shown in the picture below, paying attention to the insertion side. Then repeat the same operation for the second and the third shelf.

The side of the shelf indicated by the circle should be positioned as in the picture.

Below the final configuration is shown. Now it is possible to use the shelves to position, for example, the bags to work with TissueSAFE plus.

2.4. How to install the barcode reader

The barcode reader (fixed to its support) must be mounted on the right-hand side of the unit at the time of installation. Unscrew the two black knobs from the side panel, as shown in the image below:

Position the reader with its support and screw the knobs back on.

Finally, connect the barcode reader's cable to the USB port on the right-hand side of the TissueSAFE plus unit. The reader will be immediately ready for use.

2.5. Fume extraction system

The discharge of the TissueSAFE plus unit's vacuum pump, which is located on the back of the instrument, can be connected to a fume extraction system.

The supplied plastic tube has the following characteristics:


- Length: 3m $\pm$ 5%
- External diameter: 58 mm $\pm$ 5%
- Internal diameter: 50 mm $\pm$ 5%

The external fume extraction system must allow for an airflow of at least 12m³/h (cubic meters per hour). The unit comes with a HEPA filter and an active carbon filter, which are not accessible to the user. The first filter is designed to filter the vapours coming from the biological samples, while the second filter is designed to filter odours.

The filters are changed every 6 months during the Preventive Maintenance intervention carried out by authorized technical personnel.

Insert the discharge tube into its output on the rear panel, as shown in the following figure.

Next, connect the other end of the tube to the laboratory's fume extraction system.

2.5.1. Power supply

Power supply line needs to have a ground connection (Class I device).

Power line protection:

230V version: Miniature circuit breaker curve C, 10A, residual current device: 30mA.


Before connecting the plug to the laboratory's outlet, check to make sure that the unit's power supply voltage corresponds to that of the laboratory's network.

The use of an incorrect power supply will result in serious damage.

2.5.2. Rear panel connections

Various connections are located on the back panel of the TissueSAFE plus unit:

Description:

1. Power connector
2. Place for marking
3. Gas supply hole - not used
4. ETH/LAN connection
5. Connection for fume extraction system (for connecting the instrument to the laboratory's fume extraction system (if needed)).

The unit's label also provides important pieces of information, including the power supply voltage and the Serial Number.

Always provide the Serial Number (SN) when contacting the Technical Support Service.

Connect the power cable to the power connector shown above.


2.6. Final positioning

Once the TissueSAFE plus has been placed in its final position for use, the levers on the wheels can be lowered in order to lock them in the desired position.

Once the instrument is in its final position, lock all four wheels in place.

Before moving the instrument, make sure that all 4 of the wheels' brakes have been released.

Next, connect the power cable's plug to the laboratory's power outlet.

2.6.1. Activation/Deactivation

The main switch for activating and deactivating the instrument is located on the right panel, as shown in the figure below.

Before activating the instrument, make sure that the instructions indicated above have been carried out.

3. MILESTONE VACUUM BAGS

Milestone supplies various models of vacuum bags for use with the TissueSAFE plus unit: these bags are not supplied along with the unit, and must be ordered separately.

The TissueSAFE plus device can be used **ONLY** with the original bags supplied by Milestone s.r.l., registration certificate No.: P3H 2019/8538 for medical device: "Vacuum Bags for closing of biological specimens on TISSUE SAFE and SealSAFE".

For more details, please refer to Chapter 6 titled "Consumables".

Each bag has a white external label that can be used by the personnel to make notes.

The operator can fill out the request form and place it in the appropriate pocket.

The form can be sealed inside the pocket, thus ensuring that the accompanying document and the surgical material will be transported in a secure manner, with no possibility of errors.

4. USING THE TISSUESAFE PLUS UNIT

4.1. How to login and logout from the system


By default, the login screen appears when the device is activated, as shown in the following figure.

In order to use the main functions, the user must be registered as an administrator.


Select the "Admin" (administrator) icon, use the numeric keypad that appears to enter the password (provided on a separate document), and then press "Enter".

The main administrator screen will appear. The "Home" button can be pressed at any time from any screen in order to return to the main screen.

Temperature Data logger
NOT ENABLED

Temperature Data logger
ENABLED

Please refer to chapter 4.9 for more information on the unit's functionality with the Temperature Data Logger enabled.

The name of the active user is shown on the status bar to the left, under the date and time.

The Main Screen can be customized in the User settings by selecting Favourites and Main Functions (as explained below).

Access to the system in SERVICE mode is only permitted to authorized personnel (SERVICE mode allows for access to all the system's functions, while ADMINISTRATOR mode only allows for access to certain functions).

4.1.1. Creating a new user

Select the  Advanced Functions icon on the status bar. Next, go to the “Settings” screen and select “Users” to create a NEW USER.


The user settings screen will appear, as shown below. In order to create a new user, select “Add”. A keyboard will appear in order to enter the name of the NEW USER, which must be confirmed by pressing “Enter”. In the example below, the user name “LAB” has been entered.


Next, select and configure the parameters for the NEW USER:

- “Category”: select either “User” or “Admin” (Administrator).
- “Password”: by selecting the password box, a numeric keypad appears, which can be used to enter and confirm the password.

An example is shown in the figure.


Next, select the FUNCTIONS to be enabled for the new user:

The following boxes can be selected:


"Set Favourites as home page"	The new user will have the Favourites programs displayed on the Main Screen.
"Customize Programs"	The user will be able to create and modify programs.

By selecting the arrow shown in the figure below, it is possible to assign the users' ACCESS LEVEL upon activating the program. By default, this is set to "Ask for username and password" so that the password is always requested when the machine is turned on.

On the resulting screen shown above, the user can select the grey arrow to in order to choose whether to enable:

"Request password to execute programs"	Using a password system, the software recognizes who performs certain operations during the process (who begins or cancels a program).
--	--

The user will be asked whether they want to save the changes upon exiting the users screen.

4.2. Setting the date and time

Select the  Advanced Functions icon on the status bar. Next, go to the "Settings" screen and select "Date/time" to set the current time and date.

In the time setting window, it is possible to set the "Time Zone", thus activating the automatic update of time during the year.

After changing this setting, please reboot the system.


4.3. Owner settings

Select the  Advanced Functions icon on the status bar. Go to the “Settings” screen and select “Owner” to access the default settings.

The owner settings screen will appear, as described below.


Select “General” to modify the “Institute name” and the “Department name” where the unit is installed. An example is shown in the figure.

Select “HEPA Filter” to manage the expiration of the filter.

Below the default limits activated:

- Number of TOTAL processes to replace HEPA filter: 2000
- Number of MONTHS to replace HEPA filter: 6 months


It is strongly recommended NOT to uncheck the boxes regarding the service life of the HEPA filter.


The “Replace HEPA filter” button should be pressed once at the installation time to set zero the details about expiration date and every time the filter is replaced.

When the HEPA filter expires, and the operator presses the dedicated button to start the program, a message appears prompting the user to contact technical support. This message does not prevent the unit from being utilized, and so it is possible to run the program selected.


If you also want to active the **BLOCKING MESSAGE**, the Administrator can enable the checkboxes that by default are not activated (see below).

☒	Nr. of total processes to replace HEPA filter:	2000
☐	Nr. of total processes that stop the unit: (if HEPA filter not replaced)	0
☒	Nr. of MONTHS to replace HEPA filter:	6
☐	Nr. of MONTHS that stop the unit: (if HEPA filter not replaced)	0

Expiration date: 2018-06-22

If activated, the following blocking message appears when the operator presses the dedicated button to start the program.


When the user presses OK, the main menu is displayed and the unit can be used until the filter is replaced.

Select "Internet Connection" to modify the network settings. By selecting a numeric field, a keypad will appear in order to allow the user to enter the necessary data, which must be confirmed by pressing "Enter".

The appropriate data must be provided by the laboratory's network administrator (Primary WINS not required).

For example:

General | HEPA Filter | Internet connection | Enables | Milewatch | Print | C < >

☐ Obtain an IP address via DHCP

☒ Specify an IP address

IP address: 192 . 168 . 5 . 32


Subnet mask: 255 . 255 . 240 . 0

Default gateway: 192 . 168 . 1 . 1

Primary DNS: 192 . 168 . 1 . 6

Primary WINS: 0 . 0 . 0 . 0

Select "Enable Functions" to access the "Enable Printer" function, as shown in the figure below. The printer is enabled by default.


When "Enable Printer" is selected, at the end of each process a label containing all the relative data will be printed (for more information please refer to chapter 4.4).


If you do not want to print a label and/or paper is not available in the printer, this function must be disabled. In this manner, the software will save all the data not printed, which can be printed subsequently as soon as the printer and paper are once again available. Additional information can be found in Appendix D.1.

If the printer function is NOT enabled, the software will not display the following:


- the "Print last run" button on the main screen,

- the "Print daily report" button on the "Log book" screen, which allows the user to print the list of all the Case IDs for the programs executed on a specific day.


- The message received when a Temperature Data logger card is activated, which asks whether or not to print the activation report (as explained in chapter 4.9).


- If the laboratory uses the "Milewatch", see Operator manual MM130 to connect the unit to this software.


An example of connection is presented below.


- Select "Print" to choose the print settings. All the settings are enabled by default.

- Select "Case ID" to enable all the Case ID settings, which are explained in detail below.


- When "Enable Case ID" is selected, the software asks the user to enter a Case ID before initiating a program.

ENABLING OR DISABLING THE REQUEST FOR THE CASE ID IN ACCORDANCE WITH THE INTERNAL PROCEDURES FOR THE SPECIFIC LABORATORY (GLP - GOOD LABORATORY PRACTICES).

- When "Enable request to insert Case ID into data logger card" is selected, at the start of a program, after having entered the Case ID (if enabled), a message appears asking the user whether or not they want to enter this code on the Temperature Data Logger screen, as explained in chapter 4.9.

- When "Split case ID" is enabled, only the pertinent part of the Case ID will be taken into account. In order to do this, it is necessary to enter the desired values into the fields for the "Separator character" (character used to separate the code)


and the "Field of interest" (that which will be considered with respect to the separator character). An example is provided below:

After having entered the following values:

- 1) Separator character: ";
- 2) Field of interest: "I"

If the following Code is entered: 1111; 2222

The Case ID will become: 1111

4.4. Executing a default program

After having been cut/removed by the surgeon, the organ/sample must be placed in a Milestone vacuum bag of an appropriate size.


When executing this program following the extraction of an organ/tissue, it is recommended to remove the excess blood (and other organic fluids) with a piece of gauze. This will prevent the decay effects caused by the stagnation of these fluids.

Please refer to Appendix A for the proper use of the vacuum bags.

Follow the instructions below to execute a program.

- Once a user has logged in, the following screen appears:


- Select the "Vacuum program" icon from the main menu.

If the HEPA filter has expired, a warning message will appear, as explained in chapter 7.2.

- A "SCAN THE BAG ID" keyboard will automatically appear in order to allow the user to enter the bag ID code, as shown below.


- Once scanned, the barcode is automatically entered into the system. The bag's code is located on the lower part of the bag. The figure below shows a detailed view of the label containing the barcode: the type of bag (for example Self-standing) is indicated above, and its code in alphanumeric format is indicated below.

When the code has been read correctly, the reader emits a sound.

- If the "Enable Case ID" function (chapter 4.3) is ENABLED, when the "Bag Code" is read, another keyboard called "SCAN THE CASE ID" automatically appears. Scan the Case ID as shown below, or enter it manually using the keyboard, then press "Enter".

The maximum number of characters allowed for CASE IDs is 70.

The Case ID can contain the following characters:

- numbers: 0123456789
- uppercase letters: ABCDEFGHIJKLMNOPQRSTUVWXYZ
- lowercase letters: abcdefghijklmnopqrstuvwxyz
- symbols: !#\$%&'()*+,-.;=@ [|^_`{~
- space.


It is not possible to use characters other than those indicated.

The TissueSAFE plus keyboard does not contain: + ; = [] ` { } ~

If the Case ID reading is enabled but no code is read, the following message will appear in the Case ID field: "Not available".


If the bag codes and case ID codes are entered using the scanner (not manually), no buttons need to be pressed, and the codes are automatically included in the details of the program without having to press "Enter".

If the Temperature Data Logger is enabled and the request to automatically insert the codes into the temperature card is enabled (activation explained in chapter 4.3), when the CASE ID is read, the following question appears, which can be confirmed in order to automatically enter this code into the data logger card. For more information on how to activate the cards to record the temperature, please refer to chapter 4.9.

The code will NOT be entered if the program is interrupted by pressing "Abort" while it is in progress.

The following screen appears, which shows the Bag Code, the CASE ID (if present), and the default P1-P4 vacuum programs. An example is shown in the figure.


A message flashes above the programs indicating that the bag must be inserted before selecting a program.

- Insert the bag into the cavity. The bag containing the sample must be inserted with the opening facing the operator.

The bag must be positioned so that the lines printed on the bag itself coincide with the welding bar inside the cavity, as shown in the following figure.

Always position the uppermost line at the level of the welding bar.

Use vacuum bags of an appropriate size for the sample in question.

The use of an inappropriate bag (e.g. too small) could cause the sample to come out, and could result in errors during the vacuum and/or welding phases.


If the centreline of the sample does not coincide with the height of the welding bar, remove one or both of the spacers from the TissueSAFE plus unit's cavity.

- Select the appropriate vacuum program, from 1 to 4.

The correct program depends on the type of sample.

An “I” icon can be found below the program icons; when pressed, this icon displays the list of the types of samples for which each vacuum program can be utilized, as shown below.


In order to select the appropriate vacuum program, the user can also consult the Operator Manual's vacuum program selection guide (Appendix B).

This contains the same information that can be obtained by touching the “i” icon.

Programs P1-P4 differ by vacuum level; P1 applies the highest vacuum level, while P4 applies the lowest.

- When a vacuum program is selected (P1 for example), the following screen appears:


On the lower part of the screen, a message tells the operator to close the cover to begin the program, and to keep it pressed down for a few seconds in order to facilitate the vacuum process.

- When the cover is closed by the operator, the vacuum pump begins to create a vacuum inside the chamber.

It is important to keep the cover closed for a few seconds once the vacuum pump is activated.


The progress of the vacuum level inside the chamber is shown on the control panel under the heading “Sealing time left”.


If the "Abort" button is pressed while the vacuum process is in progress, the program is immediately stopped and the cover is automatically opened.

If the "Stop & Seal" button is pressed while the vacuum process is in progress, the program is immediately stopped, the bag is sealed, and the cover is automatically opened.

- Once the set vacuum level is reached, the bag is automatically sealed.

- Once the process has been completed, the user is notified by an audible alarm and a message stating "Run successfully completed, open cover"; the cover also opens automatically the vacuum pressure has been released.

- A label bearing the processing date is automatically printed (if the printer is enabled, please refer to chapter 4.3).


If the length of the CASE ID is greater than 14 characters, it will be automatically printed vertically on the label, otherwise it will be printed horizontally.

The data printed on the label depend on the settings entered on the Owner Data screen, as indicated in the "Print" section of chapter 4.3.


Printed data:

- Program P1-P4 (Name of the program)
- Case ID (if present)
- Bag ID
- User (User who executed the program)
- Date
- Time
- Outcome: Completed with success/Stopped & Sealed/ Aborted.

- Remove the sealed bag from the TissueSAFE plus unit.
- Tear off the label using the toothed edge.

- Remove the protective film on the back of the instrument, as shown below. The film is cut along its vertical axis to facilitate its removal. In fact, it is recommended to remove the first half first, followed by the second half.

- Attach the label to the bag.

Do not attach the label to the upper part of the bag, as it may be opened and sealed again.
 For more information on the printer, please refer to Appendix D.

If the label needs to be printed multiple times, the user can press the "Print last run" icon from the main menu.

- At this point, it is possible to note the patient's data on the bag's white label (if necessary).

- Once the surgery has been completed, the operator requesting the histological examination must fill out the "Request Form". This document must be placed in the outside pocket to the left the label where the first note is written. At this point the operator can seal the document in order to prevent it from being accidentally separated from the reference sample.

At the end of a program, it is recommended to check whether the seal is uniform along the entire surface of the bag.

If the weld is not correct:

- Open the bag with scissors.
- Clean the upper inside part of the bag.
- Repeat the vacuum program.

After having sealed the bag with a fresh sample, place it immediately in a suitable temperature-controlled refrigerator (at +4°C). An example is shown below.

The controlled-temperature system designed by Milestone s.r.l. guarantees the perfect preservation of the sample's morphology; antigenicity, DNA/RNA, and proteins are preserved in an optimal manner up to 72 hours (depending on the initial conditions and the type of sample).

The total transport time must never exceed 72 hours (including the time spent in the Operating Room or in the Histology Laboratory's refrigerator).
 In order to ensure the proper preservation of the samples' morphology, they must be stored and preserved at a temperature of +4°C±2°C.

For more details about the refrigerator models available for the TissueSAFE plus unit, please refer to chapter 5.
 The guidelines for transporting the samples can be found in Appendix C.

4.4.1. Using one bag more than once for the same patient (single use)

The software allows each bag to be sealed 3 times.


For this purpose it is extremely important for the seal to be made on the upper part of the bag, as this will allow the Pathology staff to re-utilize the bag multiple times.

The upper ends of the "Stand up" bags contain useful instructions about where to cut. For more details please refer to chapter 6.1.

The bags can be sealed three times and cut twice.

The lines to be cut are dotted.

The first seal is that which must be applied in the Operating Room, in order to leave space available for the subsequent seals.


4.4.2. Power failure during the program

In case of power failure during "Vacuum Only" program, TissueSAFE plus automatically switches off and the vacuum stays in the cavity.

When the power supply comes back on after a power failure, the system automatically switches on again, the remaining vacuum is released from the cavity and the cover can be opened.

If the power supply is not restored and / or you need to open the cover immediately, please follow the instructions below:

1. Open the cover of the bag compartment.
2. Turn the lever as shown in the picture below (red square).

3. So, the vacuum will be released and the cover can be opened.
4. Return the lever to its initial position.

Make sure the vacuum release lever is in the correct position before starting a new vacuum program.

If the power failure happened before the sealing of the bag, please follow the instructions below:


- "Vacuum Only" program: in case this type of program was performed during the power failure, repeat this program again when the power comes back on.

4.5. How to set favourite programs


This chapter contains the instructions for executing a Favourite program.

This screen can be access in two ways:


- By selecting the "Favourites" icon from the main screen,

Temperature Data logger
NOT ENABLED


Temperature Data logger
ENABLED

- After having set the Favourites screen as the main screen for a user, as explained in chapter 4.1.1.

For example, suppose that the Favourites screen has been set as the main screen for the user LAB: that which will be obtained is the following screen. The P4 default program has been set as a Favourite as an example.


When the Favourites screen has been set as the main screen, the home icon  still allows the user to return to this Screen,

while the main functions icon  allows the user to access the general screen.

From the Favourites screen, select the favourite program that you want to use following the guidelines provided in Appendix B. Once the program has been selected, follow the instructions provided in chapter 4.4 to execute the programming question.

4.6. Recovering the results of logged processes


Select the  Advanced Functions icon on the status bar, then go to the "Log book" screen.


From the log book screen, select "Saved runs".

The following screen will appear, which shows the FILES of the SAVED PROCESSES on the left, divided by date. The date is generated automatically by the software.

Select a date (29 October 2015 is used in the example below) and the processes saved on that date will appear on the right. The filename is comprised of the time in which the process was finished, followed by the letters "VO" (Vacuum Only). Click on the filename and then press "Open" to view the details (see the following figures).


The saved data are the following:

- Run name (Program name)
- User
- Start date and time
- Outcome
- Bag ID
- Case ID (if present)

Select the icon with the arrow on the lower left to return to the saved runs screen.

The user can also select the date in question and then press the "Print Daily Report" icon on the lower right: the list of the cases executed on that specific date will be printed. The list will also indicate the number of the cases that were executed without CASE CODES.

4.6.1. How to delete saved runs

Select the Date and press "Delete" to delete all the runs saved on that date (see the figure on the left).
 Select a run and press "Delete" to delete ONLY the selected file (see the figure on the right).

It is also possible to press "Delete all" to delete all the runs saved on the TissueSAFE plus unit. Press the arrow icon on the lower left to return to the main menu.

4.6.2. Exporting data from the terminal

The TissueSAFE plus unit provides all the documentation necessary for the reports on the laboratory's activities. In fact, it is possible to export the data regarding the programs executed and a document that includes the events log with all the case information.

Insert a USB flash drive into the port on the right-hand side of the instrument.

From the "Log Book" screen, select "Saved runs".

The following screen will appear. Select a data item (figure on the left) or an individual process (figure on the right) and press "Export" to export the data item or individual process in question to the USB flash drive. The user can also select "Export All" in order to export all the data present.

It is also possible to export the complete events log for the TissueSAFE plus unit.

This function is only available for Administrator and Service users.

From the "Log Book" screen, select "Events Log".

The following screen appears. Select "Export to USB". This button is only enabled when a USB flash drive is connected to the unit.

The exported data will be automatically inserted into a folder bearing the name of the TissueSAFE plus instrument followed by the instrument's serial number, as shown in the example below. This folder contains the events log "LogEvents.tissuesafeplus.s3db" and the Saved folder containing the executed programs.

The Saved folder contains the sub-folders used to separate the saved processes by execution date: the individual processes are saved within each sub-folder, and their names include the time at which they were saved.

In order to view the events log and the data regarding the executed programs on a PC after exporting them, use the LogVIEWER program provided along with the unit, and follow the instructions contained in the next chapter.


4.6.3. Installing and using LogViewer

The USB flash drive CODE 62160 supplied along with the unit contains a LogVIEWER folder. Double-click on the "SetupLogVIEWER.exe" icon to install the program.

Click "Next" on the following screens.


After having performed the steps above, the following screen will appear. Select the "Create a desktop icon" box, as shown below. Click "Next".

The "Associate s3db extension" box is selected by default, and it is recommended to maintain this setting. This will also allow the LogVIEWER program to automatically any events log files, when exported.

Finally, click "Install" to complete the installation.


This application requires Microsoft.NET Framework 4.0, which will be automatically installed on your PC if not already present. This application requires Christal Report Runtime, which will be automatically installed on your PC if not already present. At this point the LogVIEWER software is ready for use.

Insert the USB flash drive into the TissueSAFE plus unit where all the processes have been saved, and export them as described in the previous chapter.

Launch the LogViewer software (go to "Start Menu\AllPrograms\Milestone srl\Milestone LogVIEWER" or else click on the "LogVIEWER" icon on the desktop).

Insert the USB flash drive into the PC and select the icon on the upper left in the LogVIEWER window, as indicated by the green circle in the figure below.

Select the folder that was created by exporting the data from the TissueSAFE plus unit to the USB flash drive, as shown below. If you want to view the data for the exported processes, click Saved, otherwise click LogEvents to view the events log.

The Saved folder contains the sub-folders used to separate the saved processes by execution date: the individual processes are saved within each sub-folder, and their names include the time at which they were saved.

Double-click on the desired process to display the relative data. An example is shown below.

Unit name	TissueSAFE plus
Software Version	15.11
Serial number	125199002
Processing name	P1
Category	VACUUM ONLY
User	Admin
Start date	2015/12/09, at 11:01
Specimen weight	30 g
Outcome	completed with success
Bag ID	387912523026372A
Bag size	SelfStanding_LARGE

To open the Events Log, on the other hand, simply double-click on it, as shown below.

Unit name: SeaSAFE

NUMBER	TYPE	DESCRIPTION	DATE	TIME	USER	CD	LOCK	LOGN	IP
1	Aut	Reagents in	11 Nov 15	09:22:59	---	1	1	0	1409
2	Aut	User Service	11 Nov 15	09:23:14	Service (?)	1	1	0	1400
3	Aut	Access to F	11 Nov 15	09:23:27	Service (?)	1	1	0	1402
4	Aut	Access to F	11 Nov 15	09:24:20	Service (?)	1	1	0	1112
5	Aut	Set tank to F	11 Nov 15	09:24:45	Service (?)	1	1	0	1112
6	Aut	Reagents in	11 Nov 15	09:24:49	Service (?)	1	1	0	1112
7	Aut	Set tank to F	11 Nov 15	09:24:53	Service (?)	1	1	0	1114
8	Aut	Reagents in	11 Nov 15	09:24:54	Service (?)	1	1	0	1112
9	Aut	Access to F	11 Nov 15	09:25:29	Service (?)	1	1	0	0
10	Aut	Start prime	11 Nov 15	09:27:52	Service (?)	1	1	0	0
11	Aut	Start vacuum	11 Nov 15	09:37:14	Service (?)	1	1	0	0
12	Warn	Processing	11 Nov 15	09:40:15	Service (?)	1	1	0	0
13	Aut	Stop all	11 Nov 15	09:42:12	Service (?)	1	1	0	0
14	Aut	Clear all	11 Nov 15	09:43:42	Service (?)	1	1	0	0
15	Aut	Stop vacuum	11 Nov 15	09:44:25	Service (?)	1	1	0	147
16	Aut	End process	11 Nov 15	09:44:38	Service (?)	1	1	0	153
17	Aut	Add code to	11 Nov 15	09:44:39	Service (?)	1	1	0	0
18	Aut	Clear all	11 Nov 15	09:44:49	Service (?)	1	1	0	153
19	Aut	Start Temp	11 Nov 15	09:47:15	Service (?)	1	1	0	0
20	Aut	Start prime	11 Nov 15	09:47:37	Service (?)	1	1	0	0
21	Aut	Start vacuum	11 Nov 15	09:47:32	Service (?)	1	1	0	0
22	Warn	Processing	11 Nov 15	09:50:33	Service (?)	1	1	0	0
23	Aut	Stop all	11 Nov 15	09:50:33	Service (?)	1	1	0	0
24	Aut	Clear all	11 Nov 15	09:51:14	Service (?)	1	1	0	0
25	Aut	Stop vacuum	11 Nov 15	09:52:27	Service (?)	1	1	0	152
26	Aut	End process	11 Nov 15	09:52:40	Service (?)	1	1	0	15


4.7. How to modify a program


This function does NOT create a NEW PROGRAM. It is necessary to select the standard program that best meets your needs, make a copy of it, and modify the details.

Before modifying a program, contact the application apartment at the following email address: application@milestonemedsl.com.

Select the  Advanced Functions icon on the status bar, then go to the "Programs" screen.

The following screen will appear with the current list of programs, including the pre-set protocols and the protocols created by the users (if present).


In order to modify a default protocol, select it, click “Copy”, rename the program, and modify its details.

Modification of a user's protocol does NOT require copying it before making changes required.

When “Copy” is selected, a program called “Copy of...” is automatically created: clicking on this name will call up a keyboard that will allow you to enter the desired name.

In order to modify the vacuum level, select “Modify vacuum level”: the program's vacuum value will be displayed. Enter the desired vacuum level and click on the arrow to save the new program.

It is not possible to create a program with a vacuum level greater than 600 mbar, otherwise it wouldn't be possible to seal the bag. Such a value would nevertheless correspond to a seal “without a vacuum” inside the bag.

When modifying a protocol created by a user there is no need to make a copy beforehand.

The “Modify” and “Delete” buttons are only available for programs that have already been copied. For security purposes, the protocols pre-set by the manufacturer cannot be modified or deleted.

The “Modify” function can be used to modify the vacuum details of a program previously created by the user, but not its name.


4.8. How to set favourite programs

The FAVOURITE programs can be included on the MAIN SCREEN in the following manner:


- Access the advanced functions by clicking on the relative icon on the status bar. Next, go to the “Settings” screen and select “Set favourites”.

- Click on the "Not Defined" button under the "Choose the Favourite to set" heading.

- The list of current programs will appear. Select the desired program. The colour of the background can be changed by selecting "Choose the button colour", as shown in the following image, for example.


- Click "Save favourite" to save the program.
- The program's icon will now be shown on the favourites screen. The favourites screen can be set as the main screen, as explained in chapter 4.1.1.


4.9. Temperature Monitoring System

When the TissueSAFE plus instrument is installed in the Operating Unit and is equipped with an integrated card reader, it can be used to activate the cards in order to monitor the transport temperature of the samples to Pathology.

In this case, the main screen of the TissueSAFE plus unit's software contains the "Temperature logger" icon in order to access this function.

The following image shows the screen that is opened by clicking on the "Temperature logger" icon.

The upper portion of the screen contains the following information:

- Unit ID: indicates the instrument's serial number (SN)
- Status: indicates the status of the card to be activated when it is inserted into the TissueSAFE plus unit
- Battery level: indicates the level of the card's battery, which must be "Good" in order to be utilized
- Firmware: indicates the card's firmware version.


When the battery level is "Low", the user is advised not to use the card as the remaining power may not be sufficient to complete the temperature logging operations.

Prior to activating the card, it is possible to enter the codes of the samples that will be sent along with the card itself.

This function is not mandatory.


The codes can be entered in two manners:

1. Manually, by clicking on the "Modify bar codes" icon . The following screen will appear.


The codes can be automatically inserted on this screen by scanning them with the barcode reader.


Otherwise, they can be inserted manually by clicking on "Add code" and using the keyboard that appears. The final result is the following (for example):

On the same screen it is also possible to delete any codes entered by mistake using the "Delete code" and "Clear list" buttons. Click on the arrow on the lower left to return to the temperature logger's main screen.

2. Automatically, by activating the request to automatically insert the codes when executing a vacuum program, as explained in the "Case ID" section of chapter 4.3.
If this request is activated, after entering the Case ID at the start of a vacuum program, a message will appear asking the user whether or not they want to automatically enter this code on the card for temperature logging.

It will therefore no longer be necessary to enter the "Modify barcodes" section as explained under point 1, because in automatic mode the codes will already be entered and it will only be necessary to activate the card for the temperature logging to be carried out. This section can nevertheless be accessed if they need to be deleted.

If the vacuum program is interrupted by clicking "Abort", the code will not be entered even if the request to enter the codes has been confirmed.

It is not necessary to keep the card inserted into the TissueSAFE plus unit during the automatic insertion of the codes, as the codes will be saved on the dedicated screen (it is nevertheless important not to shut off the instrument in order to avoid losing the codes). For example, the card could be kept in the refrigerator and only inserted into the TissueSAFE plus unit at the time of activation.

Before activating the card, it is necessary to select the time interval that must pass between one temperature measurement and the next.

This value must be selected based on the duration of the logging operations, and therefore depends on the distance between the Operating Unit and the Pathology Unit. The possible values are the following:

- 2 minutes: up to 30 hours of logging
- 5 minutes: up to 3 days of logging
- 10 minutes: up to 6 days of logging

Read temperature every:


- ☐ 2 min (up to 30 hours)
- ☒ 5 min (up to 3 days)
- ☐ 10 min (up to 6 days)

Activate the card by clicking on the "Start" button.

Once the card has been activated, a sound is emitted to confirm its activation, and the card's status changes to "Running".

Unit ID: TissueSAFE plus_150399001
 Status: Running, started on Tue, 3-Nov-2015, at 11:32
 Battery level:  Good Card firmware: 4

If the printer is enabled, a message will appear asking whether the user wants to print the card's activation report.

Printed data:

- SN (the unit's Serial Number)
- Date (Activation Date)
- Time (Activation Time)
- CASE ID list (the list of case IDs included on the card) – only if present.

If necessary, the card can be activated with a delay by clicking the **Set start date** button.

The panel shown in the following figure on the right appears. Set the desired start date and time and click “Start”.

The status of the activated card will change to “Delayed to...”.

Unit ID: TissueSAFE plus_150399001
Status: Delayed start on Wed, 4-Nov-2015, at 11:51
Battery level:  Good **Card firmware:** 4

4.10. Sample transportation

The fresh vacuum-packed samples must be stored at +4°C.

Milestone s.r.l. offers a temperature monitoring system for keeping the temperature under control, in accordance with the sample transportation guidelines. **The sample transportation guidelines are given in Appendix C.**

Once the courier is ready to transport the samples to the Pathology unit, the vacuum-packed samples must be placed in the refrigerator box with 1, 2, 3, or 4 cool packs (based on the space available inside) at +4°C, after which the box must be closed and transported (figures below).

In order to ensure the optimal preservation of the samples, it is recommended to keep the cooling blocks in a refrigerator near the TissueSAFE plus unit, at a temperature of +4°C ± 2°C.

Do not place or transport cool packs at temperatures lower than +2°C inside the portable cooler. Frozen packs can damage the samples during transportation.

The cool packs are extremely important for maintaining the temperature inside the refrigerator box. Never transport any samples in the box without cool packs present.

Once the transport operations have been completed, the samples are ready to be cut open.

4.11. How to use the bag sealer

The bag sealer can be used to keep closed temporary a bag containing or not a specimen, dirty or not by fixative.

How to use the bag sealer:


- Bend the upper part of the bag as shown into the picture.

- Insert the yellow spindle of the bag sealer into the bended part of the bag.

- Then push the bag sealer along the bended part of the bag.

- The bag sealer will keep closed all the opened side of the bag.

Below the final result.

It is possible to use the handle of the bag sealer to move a bag inside the laboratory but it is required to keep one hand below the base of the bag for safety.

5. COMPONENTS AND ACCESSORIES

Table 2 provides the characteristics of components and accessories which are included in TissueSAFE plus:

Table 2

Technical specification (TS)	TS value / parameter
Components	
Discharge tube	Intended use - Thermoplastic rubber hose with coated steel spiral is used for directing the exhaust vapor flow Q-ty 1 pc. Internal diameter = 50mm$\pm$5% External diameter = 58mm$\pm$5% Weight = 340g per meter$\pm$10% Temperature range: from -25° C to +125° C Length - 3 m$\pm$10% 
Metal strap	Dimensions = 60x80x13 mm$\pm$5% Weight: 31 g $\pm$10% Intended use - Stainless steel hose strap to secure the clamp Q-ty 2 pcs 

Table 2 (continuation)

Technical specification	TS value / parameter																								
Main power cable	Dimensions (length x diameter) : 1.8 m (±5%) x 4.8 mm (±5%) Weight: 205 (±10) g Intended use - Power cable Schuko-VDE for powering the device Q-ty 1 pc.																								
																									
Bar code reader mounted on its support	Q-ty 1 pc. Color - Black Dimensions: 181 x 100 x 71 mm (±5%) Weight: 195.6 g ±10% Input voltage: 4.2-5.25 V DC Operating rated currents: 160 mA Max. current consumed: up to 350 mA Current consumed in idle mode: 65 ma Performance characteristics: - Light Source: LEDs - Max reader roll angle: up to±180°; - Max angle of skew/rotation of the reader relative to its support: ± 40°; - Print contrast minimum: 25% minimum reflectance; - Drop specification: Scanner withstands 18 drops from 1.8 meters to concrete; - Ambient Light Immunity: up to 100,000 Lux - Contaminants: Spray/rain dust/particulates: IEC 529-IP52; - ASD Level: 16 kV; - Depth of field (see the Table below):																								
	Laser's characteristics: - Power consumption: up to 1 mWt, - Max. wave length (λ) = up to 650 nm - Pulse width = 12 mcs - Laser device of II class - Beam divergence 8,4° x 8,1° ("plus")																								
	<table><tr><th colspan="3">Depth of Field (Typical) *</th></tr><tr><th>Symbol</th><th>SR:</th><th>HD:</th></tr><tr><td>Code 39</td><td>5 mil: 1.6"-7.5" (4,0-19 cm) 10 mil: 0.4" - 11.8" (1,0-30 cm) 20 mil: up to 17.7" (up to 45 cm)</td><td>3 mil: 0.9"-3.6" (2,4-9,1 cm) 5 mil: 0.3"-4.5" (0,8-11,3 cm)</td></tr><tr><td>EAN</td><td>7,5 mil: 0.5"-10.6" (2,0-27 cm) 13 mil: 0.6"-15.7" (1,5-40 cm)</td><td>7,5 mil: 0.5" (0-12,7 cm) 13 mil: 43"-6.8" (1,1-17,2 cm)</td></tr><tr><td>PDF-417</td><td>6,6 mil: 1.0"-5.9" (2,5-15 cm) 10 mil: 0.2"-8.6" (0,5-22 cm) 15 mil: 0.6"-13.4" (1,5-34 cm)</td><td>4 mil: 0.7"-2.7" (1,8-6,8 cm) 6,6 mil: 0.1"-4.4" (0,1-11,2 cm) 10 mil: 0"-5.6" (0-14,3 cm)</td></tr><tr><td>DataMatrix</td><td>10 mil: 0.8" to 6.3" (2,0-16 cm) 15 mil: 0" to 9.3" (0-23,6 cm)</td><td>5 mil: 1.1"-2.4" (2,8-6,1 cm)</td></tr><tr><td>QR-code</td><td>10 mil: 1.2" to 4.9" (3,0-12,5 cm) 15 mil: 0.4 to 7.5" (1,0-19 cm)</td><td>6,7 mil: 0.8"-1.7" (2,1-4,2 cm)</td></tr><tr><td>Minimum element width</td><td>Standard range: 1D Min Resolution: 4 mil PDF-417 Min Resolution: 5 mil DataMatrix Min Resolution: 7 mil</td><td>High Density: 1D Min Resolution: 2,5 mil PDF-417 Min Resolution: 4 mil DataMatrix Min Resolution: 5 mil</td></tr></table>		Depth of Field (Typical) *			Symbol	SR:	HD:	Code 39	5 mil: 1.6"-7.5" (4,0-19 cm) 10 mil: 0.4" - 11.8" (1,0-30 cm) 20 mil: up to 17.7" (up to 45 cm)	3 mil: 0.9"-3.6" (2,4-9,1 cm) 5 mil: 0.3"-4.5" (0,8-11,3 cm)	EAN	7,5 mil: 0.5"-10.6" (2,0-27 cm) 13 mil: 0.6"-15.7" (1,5-40 cm)	7,5 mil: 0.5" (0-12,7 cm) 13 mil: 43"-6.8" (1,1-17,2 cm)	PDF-417	6,6 mil: 1.0"-5.9" (2,5-15 cm) 10 mil: 0.2"-8.6" (0,5-22 cm) 15 mil: 0.6"-13.4" (1,5-34 cm)	4 mil: 0.7"-2.7" (1,8-6,8 cm) 6,6 mil: 0.1"-4.4" (0,1-11,2 cm) 10 mil: 0"-5.6" (0-14,3 cm)	DataMatrix	10 mil: 0.8" to 6.3" (2,0-16 cm) 15 mil: 0" to 9.3" (0-23,6 cm)	5 mil: 1.1"-2.4" (2,8-6,1 cm)	QR-code	10 mil: 1.2" to 4.9" (3,0-12,5 cm) 15 mil: 0.4 to 7.5" (1,0-19 cm)	6,7 mil: 0.8"-1.7" (2,1-4,2 cm)	Minimum element width	Standard range: 1D Min Resolution: 4 mil PDF-417 Min Resolution: 5 mil DataMatrix Min Resolution: 7 mil
Depth of Field (Typical) *																									
Symbol	SR:	HD:																							
Code 39	5 mil: 1.6"-7.5" (4,0-19 cm) 10 mil: 0.4" - 11.8" (1,0-30 cm) 20 mil: up to 17.7" (up to 45 cm)	3 mil: 0.9"-3.6" (2,4-9,1 cm) 5 mil: 0.3"-4.5" (0,8-11,3 cm)																							
EAN	7,5 mil: 0.5"-10.6" (2,0-27 cm) 13 mil: 0.6"-15.7" (1,5-40 cm)	7,5 mil: 0.5" (0-12,7 cm) 13 mil: 43"-6.8" (1,1-17,2 cm)																							
PDF-417	6,6 mil: 1.0"-5.9" (2,5-15 cm) 10 mil: 0.2"-8.6" (0,5-22 cm) 15 mil: 0.6"-13.4" (1,5-34 cm)	4 mil: 0.7"-2.7" (1,8-6,8 cm) 6,6 mil: 0.1"-4.4" (0,1-11,2 cm) 10 mil: 0"-5.6" (0-14,3 cm)																							
DataMatrix	10 mil: 0.8" to 6.3" (2,0-16 cm) 15 mil: 0" to 9.3" (0-23,6 cm)	5 mil: 1.1"-2.4" (2,8-6,1 cm)																							
QR-code	10 mil: 1.2" to 4.9" (3,0-12,5 cm) 15 mil: 0.4 to 7.5" (1,0-19 cm)	6,7 mil: 0.8"-1.7" (2,1-4,2 cm)																							
Minimum element width	Standard range: 1D Min Resolution: 4 mil PDF-417 Min Resolution: 5 mil DataMatrix Min Resolution: 7 mil	High Density: 1D Min Resolution: 2,5 mil PDF-417 Min Resolution: 4 mil DataMatrix Min Resolution: 5 mil																							
																									

Table 2 (continuation)

Technical specification	TS value / parameter
Plastic spacers positioned in the cavity	Dimensions: 350x340 mm ($\pm 5\%$) Weight: 3 kg $\pm 5\%$ Purpose - Holding the bag inside the chamber Q-ty 2 pcs 
Shelves to be inserted into the holder for bags	Dimensions: 412x335 mm ($\pm 5\%$) Weight: 1.94 kg $\pm 5\%$ Intended use - For storage of bags Q-ty 3 pcs 
USB flash drive containing "LogVIEWER" software	Dimensions: Length - 66 mm $\pm 5\%$, Width - 20 mm $\pm 5\%$ Weight, 5 g $\pm 10\%$ Intended use - USB flash drive (8 GB) containing «LogVIEWER» software USB flash drive also contains the Operator Manual in electronic format. Q-ty 1 pc. SW version 20.07 Release date: 13.07.2020 LogVIEWER functionality: - unit control; - optimization of user interaction with the device; - display of tasks, parameters and results of the device operation. - reading reports on the TissueSAFE plus operation. The reports are downloaded separately from TissueSAFE plus by a specialist using a USB key. The report includes the event log and saved processes. 

Table 2 (continuation)

Technical specification	TS value / parameter
Rolls of printer paper	Paper type: adhesive Dimensions: Outer diameter 55 mm $\pm$5% / Internal diameter 18 mm $\pm$5% / Height 57 mm $\pm$5% Weight: 11 g $\pm$10% Intended use - for printing the labels to vacuum bags Q-ty 10 pcs Storage conditions (without packaging): at +10°C... +40°C. Do not expose to sunlight 
Stylus for working with the touch screen terminal	Dimensions length: 150 mm $\pm$5% / diameter 7.5 mm $\pm$5% Weight: 11 g $\pm$10% Intended use - for working with the touch screen terminal Q-ty 1 pc. 
Bag sealer with handle	Dimensions: 350x65x11 mm $\pm$5% Weight: 50 g $\pm$10% Intended use - For fixation of bags Q-ty 12 pcs 

Table 2 (continuation)

Technical specification	TS value / parameter
Kit Refrigerator box in PS, 32 liters complete with N.4 cool packs (1 pc)	Intended use - samples packed in vacuum bags are placed in a refrigerator box. Cool packs and a temperature monitoring system (if present) are also placed inside the refrigerator box. If the temperature regime inside the refrigerator box is maintained at $+2^{\circ}\text{C} \dots +4^{\circ}\text{C}$, the samples are preserved their biological properties for 72 hours. Samples in the refrigerator box can be transported or stored for 72 hours. Each cool pack is filled with 250 g ethylene glycol (4 pcs). Warning! After sealing, biological samples should be stored at $+2^{\circ}\text{C} \dots +4^{\circ}\text{C}$. For this purpose, the Kit Refrigerator box is included in the delivery set of the instrument. Without the Kit Refrigerator box, the function of the instrument, transportation and/or storage of samples is not implementable. This Kit Refrigerator box is delivered and used only with TissueSAFE plus, independent use of the box is unacceptable. Information on the devices for carrying (transportation): The Kit Refrigerator box is suitable for transportation and storage in accordance with cl. 10 of this Manual. Information on ways to maintain and control the temperature range: The temperature regime inside the refrigerator box is maintained by means of 4 cool packs. To control of the temperature range in the Kit Refrigerator box can be used the data logger cards, which is a part of Temperature monitoring system. One data logger card together with the packed biological specimens is placed in the center of the container and records the temperature and its changes during the transportation and/or storage of specimens. Information from the data logger card is read by using of the data logger card reader and visualized as a temperature graph. If the temperature during transportation and / or storage was outside the normal range ($+2^{\circ}\text{C} \dots +4^{\circ}\text{C}$) the samples are not suitable for further research. And temperature indicators WarmMark by ShockWatch, Inc., USA (registration certificate No. $\Phi\text{C}3$ 2008/02988, version: WarmMark Shot Run, activation temperature: $+5^{\circ}\text{C}$) can be used to control the temperature. Information on design: Kit Refrigerator box in PS, 32 liters, rectangular shape, with a cover, it has a handle for ease of transfer. Dimensions: $\text{WxDxH } 490 \times 310 \times 420 \text{ mm} \pm 5\%$ Weight: $7 \text{ kg} \pm 10\%$ Q-ty 1 pc Cool packs of rectangular shape Cool pack dimensions: $21 \times 18 \times 2.5 \text{ cm} \pm 5\%$ Weight: $0.4 \text{ kg} \pm 10\%$ Q-ty: 4 pcs Information on the material used for manufacturing: Material of the Refrigerator box: polypropylene (PP), Material of the Refrigerator box's cover: polystyrene (PS), Insulating material of the Refrigerator box: polyurethane (PUR), Material of the cool packs - high-density polyethylene (HDPE), filling material - ethylene glycol Principle of operation and placement of cool packs and biological samples (quantity, location in the box) to achieve the specified temperature range "$+2^{\circ}\text{C} \dots +4^{\circ}\text{C}$ for 72 hours": biological samples sealed under vacuum must be placed in the refrigerator box together with 4 cool packs at $+4^{\circ}\text{C}$, close tightly. Cool packs must be placed on each side of the box. The quantity of sealed bags with biological samples that can be placed in the refrigerator box depends on the sample size. It is recommended to fill the box no more than 2/3 of the volume. For optimal archiving of the samples Milestone s.r.l. recommends storing the cool packs in a refrigerator next to TissueSAFE plus at $+4^{\circ}\text{C} \pm 2^{\circ}\text{C}$. Do not store or transport cool packs in the refrigerator box at temperatures below $+2^{\circ}\text{C}$. Frozen cool packs can damage biological samples during transportation. The total transportation time should not exceed 72 hours (including time spent in the operating room or the refrigerator of the histological laboratory). Frequency of use details: Reusable Information on preparation of the refrigerator box and cool packs for their subsequent use:

	1. Pre-cool the cool packs in the refrigerator to +2°C...+4°C. Approximate cooling time to the desired temperature is 8-12 hours. Place 4 cool packs in the refrigerator box with sealed samples. 2. Disinfect the box and cool packs with disinfectants specified in section 9. It is necessary to clean them after each use. Operating conditions of the refrigerator box: temperature range: + 2 °C.. + 30 °C, and relative humidity: 10–80%, non-condensing. Manufacturer specific requirements. Classification of the refrigerator box depending on the perceived mechanical impacts: Group 3 according to GOST R 50444 Parameters of resistance of the refrigerator box to the specified impacts: The refrigerator box is enough resistant to shocks and vibrations. Resistant to mechanical impacts during transportation, to vibrations and impacts. It is resistant to exposure of biological fluids and secretions of body tissues. It is resistant to polyphenolic and alcohol-based disinfectants and cleansers.
Temperature monitoring system	Purpose: for temperature reading and monitoring Delivery set: temperature control kit composed of coding device, data logger card reader and kit of 10 data logger cards Main characteristics: • Increment: $\pm 1^{\circ}\text{C}$ • Temperature is read each 2/5/10 minutes • Operating temperature: from -10°C to 40°C • Service life of the battery: 2 years. • Dimensions of card reader (USB cable of 1.5 m length): $115 \times 45 \times 27 \text{ mm} \pm 5\%$ • Dimensions of data logger card: $80 \times 40 \times 12 \text{ mm} \pm 5\%$ • Maximum number of characters in case ID: up to 240 
Accessories	
HEPA filter	Type/Model: 68290 Dimensions: length $14.7 \text{ cm} \pm 5\%$ x diameter $5.5 \text{ cm} \pm 5\%$, weight $132 \text{ g} \pm 10\%$ Q-ty 1 pc Intended use: for filtering vapors from biological samples. 

Carbon filter	Type/Model: 68274 Dimensions: 15 cm x 15 cm x 6.8 cm $\pm 5\%$, weight 1300 g $\pm 10\%$ Q-ty 1 pc Intended use: for filtering odors. 
---------------	---

6. ADDITIONAL DEVICES THAT ARE NOT PART OF THE INSTRUMENT, BUT PROVIDED FOR CO-USE

6.1. Milestone Vacuum Bags

Combined with the dedicated welding bar, the particular design and special multilayer construction of Milestone vacuum bags provide for a perfect double seal that even holds in the case of accidental falls or external compression. Milestone vacuum bags must therefore only be utilized with the TissueSAFE plus instrument. The use of these bags in instruments other than that which is indicated, or the use of bags not supplied by Milestone in the instruments indicated above, could compromise their safety, thus resulting in biological risks for the operator and the risk of compromising the quality of the samples themselves, while at the same time voiding the manufacturer's warranty and the unit's CE-IVD certification.

Every bag is marked with the "Single use only" symbol (see the circle in the figure below), which indicates that the bag may only be used for the same sample.

Milestone s.r.l. offers six types of vacuum bags, two of which are classified as stand up.

All the types of vacuum bags offered by Milestone are compatible with the TissueSAFE plus instrument.

Information on how to use the vacuum bags on TissueSAFE plus is provided in section 4.4 and Appendix A.

6.1.1. Storing the vacuum bags

Store in a cool dry place, away from any heat sources or direct sunlight.

6.1.2. Stand up vacuum bags

Milestone offers two different sizes of self-standing vacuum bags, which only differ in terms of dimensions. An image is provided below.


Vacuum bag SU large (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 330mm x 440mm (± 7 mm) Q-ty 200 pcs	Purpose: for primary packaging of biological samples for further vacuum sealing in TissueSAFE plus
Vacuum bag SU small (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 220mm x 310mm (± 7 mm) Q-ty 200 pcs	Registration certificate No. P3H 2019/8538 for Vacuum Bags for closing of biological specimens on TISSUE SAFE and SealSAFE.

Characteristics:

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient)
- Due to the special design with the stand up bottom, the vacuum bag is stand up and can stand alone on a table or other horizontal surface
- High mechanical strength
- With applied sealable and re-openable document pocket
- Solid lines indicating where to seal
- With dotted lines indicating where to cut
- Double barrier layer of material BOPA/PE (polyamide and polyethylene)
- With indication of the batch number and expiry date
- Enveloped in packs of 50 pcs

6.1.3. Vacuum bags (not stand up)

Available in four models:

Vacuum bag small (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 250mm x 300mm (± 7 mm) Q-ty 200 pcs	Purpose: for primary packaging of biological samples for further vacuum sealing in TissueSAFE plus Registration certificate No. P3H 2019/8538 for Vacuum Bags for closing of biological specimens on TISSUE SAFE and SealSAFE.
Vacuum bag medium (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 330mm x 385mm (± 7 mm) Q-ty 200 pcs	
Vacuum bag large (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 330mm x 555mm (± 7 mm) Q-ty 200 pcs	
Vacuum bag extra large (200 pcs)	Thickness of the vacuum bag side 145 μm ($\pm 10\%$), External dimensions: 330mm x 555mm (± 7 mm) Q-ty 200 pcs	

The four types of vacuum bags described above differ only in size. Their common characteristics:

- Disposable bags (if necessary prescribed the sealing of the bag up to 4 times (maximum) for the same specimen from the same patient)
- High mechanical strength
- With applied sealable and re-openable document pocket
- Solid lines indicating where to seal
- Double barrier layer of material BOPA/PE (polyamide and polyethylene)
- With indication of the batch number and expiry date

6.2. Temperature indicators

In addition to the temperature data logger cards included in the unit delivery set for reading and controlling the temperature, the following temperature indicators – WarmMark by ShockWatch, Inc., USA, registration certificate No. $\Phi\text{C}3$ 2008/02988 for a medical device: “Temperature indicators: WarmMark, ColdMark” - can be used to control the temperature regime in a kit refrigerator box.

Version: WarmMark Shot Run

Activation temperature: $+5^{\circ}\text{C}$

Dimensions: 46x19x2.5 mm ($\pm 5\%$)

Usage of indicators.

Installation of WarmMark indicators: before installing the indicator, it is necessary to keep it for an hour at temperature at least 5°C below the activation temperature for complete solidification of the dye. Then, without allowing the cooled indicator to heat up, remove the base surface protecting the adhesive layer and press the adhesive layer against a clean, dry surface, or place it simply next to the controlled item. Bringing the indicator into working condition: break the label along the perforated line and

pull it to remove the protective film. You can also bring the indicator into working condition before gluing, but it should not have time to heat up before installation.

7. MAINTENANCE

The maintenance of the TissueSAFE plus unit must be performed on two levels.

- 1) The first level is the laboratory's responsibility, and consists of simple maintenance and routine cleaning operations, which are extremely useful for preventing problems.
- 2) The second level is more technically demanding and requires specialized skills typically not required of laboratory personnel.

Attention: please contact your local representative or the Milestone Customer Support Service directly (customersupport@milestonemedsrl.com) for more detailed information on how to perform Routine Maintenance and regular inspections in order to keep your TissueSAFE plus unit in perfect working order.

7.1. Operating room personnel

The users are required to perform minimal maintenance activities in order to the TissueSAFE plus unit's optimal performance. The following operations are recommended:

7.1.1. Daily

Decontaminate and clean the instrument according to the procedure described in chapter 9.

Make sure that the TissueSAFE plus unit has been cleaned and the contaminated; do not use abrasive substances or sharp objects that could damage the casing.

Check the cover and edges of the vacuum compartment.

Regularly remove any traces of organic materials in order to ensure the cover's perfect closure.


The cover's imperfect closure could compromise the unit's ability to generate a vacuum.

7.1.2. Weekly

The vacuum pump must be cleaned at least twice a week. No vacuum bags are required for this operation.


The procedure only takes 10 minutes and will increase the service life of the vacuum pump.

From the main menu, click on the "Pump Cleaning" icon.

A message will appear prompting the user to close the cover and hold it down for a few seconds in order to initiate the procedure. Once the cover is closed, the process begins automatically, with a vacuum being generated in order to thoroughly clean the pump. While the program is in progress (10 minutes), the screen displays a progress bar to show the time remaining.

The procedure can be cancelled at any time by pressing "Abort".

Once complete, the unit automatically returns to the main menu.

7.1.3. As required

When the paper runs out, replace the roll as explained in appendix D.1.2.

7.1.4. Yearly

Contact your technical support for a preventative maintenance contract.

7.2. Service staff

Please contact your local representative or the Milestone Customer Support Service directly (customersupport@milestonemedsrl.com) for more detailed information on how to perform Routine Maintenance and regular inspections in order to keep your TissueSAFE plus unit in perfect working order.


The unit comes with a HEPA filter, which are not accessible to the user, and it is designed to filter the vapours coming from the biological samples.

The filters are changed every 6 months during the Preventive Maintenance intervention carried out by authorized and trained technical personnel.

The HEPA filter is also managed by the software.


When the HEPA filter expires, a message appears when the program is launched prompting the user to contact technical support. This message does not prevent the unit from being utilized.

By default, this message appears once the instrument has executed 2000 cycles.

If the control function is enabled, a message will appear after a certain number of processes are executed with the filter expired, notifying the user that the unit will no longer be able to be utilized until the filters are replaced.

By default, this control function is not enabled.

Repair of the device is carried out only by Milestone S.r.l. It is strictly forbidden to open the device without an authorized representative of the manufacturer. The manufacturer warrants the continued safety of the equipment after maintenance and checks. The fuses integrated in the device are overcurrent fuses, and they are not subject to replacement by users. If you need to replace the fuse, please notify the Milestone service engineer.

8. WASTE DISPOSAL

8.1. Disposal of vacuum bags

The Vacuum bags supplied by Milestone s.r.l. should be treated as biohazard material after their use.

Waste material should be disposed of in an approved incinerator in compliance with all federal, provincial and local government regulations.

For the Russian Federation it is recommended to be disposed of vacuum bags in accordance with the rules implemented by "SanPiN 2.1.3684-21 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class "B" (epidemiologically hazardous waste).

8.2. Waste of the unit, its components and accessories

Recommendations for removal of the equipment for disposal:

- 1) Shut off the unit. Disconnect the unit from the electrical outlet.
- 2) Disconnect all the components and accessories.
- 3) Clean and disinfect the unit, its components and accessories in accordance with Chapter 9 of this manual.
- 4) Eliminate all of the equipment's electrical and electronic components.

Disposal and environmental requirements

Disposal of the unit, its components and accessories must be carried out in accordance with the rules of collection, recording and disposal established by the authorized federal executive body and provided for electronic devices.

It is not allowed to dispose it of as household waste.

Before disposal, the unit and its components are subject to mandatory decontamination (disinfection) / neutralization. The choice of the method of decontamination / neutralization is determined by the capabilities of the organization carrying out the medical activities, and it is implemented when developing a scheme for handling of medical waste.

The unit is subject to disposal in the following cases:

- end of service life;
- in cases endangering the life and health of the service personnel or the patient.

All packaging, including transport packaging, must be also disposed of.

For the Russian Federation it is recommended to be disposed of the unit in accordance with the rules implemented by "SanPiN 2.1.3684-21 (Sanitary and epidemiological requirements for the medical waste handling)" for the medical waste referred to class "A" (non-hazardous waste), excluding electronic parts. Electronic parts must be disposed of in accordance with local specified regulations with respect to household and industrial waste (in particular, office and computer equipment).

9. DECONTAMINATION AND CLEANING

The procedure described in this chapter must be performed every day.

For safety reasons, any operations that do not involve the use of personal protection equipment must be avoided.

The instructions reiterate the user's responsibility to properly decontaminate the unit in the event of any spills of hazardous materials on or within the unit itself.

The instructions reiterate the user's responsibility NOT to use decontamination or cleaning agents that could pose a hazard by reacting with parts of the instrument itself or any contaminated materials therein.

This instrument is an in-vitro medical device and it is usually installed in a laboratory where specimens and other biological tissues are present. For your safety, it is therefore required to clean and disinfect before entering in contact with it. You must wear gloves when operating the system.

Do not remove the side panels during cleaning procedures.

Disconnect the device from its electrical power supply before performing any operations.
Do not use direct or high-pressure water to clean the instrument.

For safety reasons, any operations that do not involve the use of personal protection equipment must be avoided.

For more information on TissueSAFE plus disinfection and cleaning see chapter 9.2

9.1. Safety precautions and regulations

Any personnel performing or otherwise involved in the decontamination of the TissueSAFE plus unit, its components and accessories must make use of suitable personal protective equipment, namely:

- Single-use HI-RISK type protective gloves (100% HI-RISK latex or equivalent material). Personal protection device compliant with the EN 374 standard and the 89/686/EEC Directive.

- An antibacterial mask for surgical/operating room use and a plastic apron. Personal protection device compliant with the EN 149 standard Grade FFP3D and the 89/686/EEC Directive.

- Protective eyewear. Personal protection device compliant with the EN 166 standard and the 89/686/EEC Directive.

9.2. Procedure for decontaminating and cleaning the TissueSAFE plus unit

9.2.1. Detergents

The following instructions must be followed scrupulously in order to ensure the proper decontamination and cleaning of the unit.

Do not remove the unit's side panels during cleaning operations.

Do not use any procedures other than those described in this chapter for the decontamination and cleaning operations.

Disconnect the unit from its electrical power supply before performing any operations. Do not use direct or high-pressure water to clean the instrument.

It is recommended to use polyphenol-based disinfectants for the decontamination operations and alcohol-based products for the cleaning operations (except the cover, as explained below).

DO NOT use alcohol-based detergents to clean the cover.


If you have questions or for more information, please contact our specialists at the following email address: application@milestonemedsrl.com.

In accordance with the rules established by MU 287-113 "Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical devices", it is allowed to clean the working surfaces with a 3% solution of hydrogen peroxide (for example, produced by JSC "Samamedprom", Russia) with the addition of 0,5% "Lotos" type detergent (for example, produced by "RHP" LLC, Russia). It is allowed to use cleaning and disinfecting agents from other manufacturers that are similar in composition and allowed for circulation in the territory of the Russian Federation.

9.2.2. Procedure for decontaminating and cleaning the unit

1) Spacers

- Extract the spacers from the work area.
- Spray with a disinfectant solution and let stand.
- Rinse with running water to remove the disinfectant.
- Wash with an appropriate detergent and, if necessary, remove any traces of organic materials with a non-abrasive instrument.
- Rinse with running water to remove the detergent.
- Dry with a clean cloth or paper towels.

2) Cover

For cleaning operations it is recommended to use water and soap or any liquid detergent. Do not use alcohol or acids.

- Spray a disinfectant solution onto the plexiglass cover (both inside and outside) and let stand.
- Clean the surface with a water-soaked cloth to remove the disinfectant.
- Use an appropriate detergent to clean the cover, and remove any organic residues that may be present a dry cloth or paper towel.
- Clean the surface with a water-soaked cloth to remove the detergent.
- Dry with a clean cloth or paper towels.
- Use a dry cloth or paper towels to thoroughly wipe off the red silicone closure bar (indicated with the letter A in the figure below).
- Execute the default P1 program without using a bag in order to thoroughly dry the silicone bar. Run the program a second time if necessary, and then wipe it off again with paper towels.
- After having run the vacuum program, dry the welding bar in the cavity with paper towels as well.

Throughout the procedure described above for the cover, pay particular attention to the black rubber seal on the inner surface (marked with B on the previous photo). It must be perfectly clean and dry in order to properly adhere to the steel surface of the operating compartment during the execution of a vacuum program.

3) Cavity

- Spray with a disinfectant solution and let stand.
- Clean the surface with a water-soaked cloth to remove the disinfectant.
- Use an appropriate detergent to clean the cavity, and remove any organic residues that may be present using a dry cloth or paper towel.
- Clean the surface with a water-soaked cloth to remove the detergent.
- Dry the surface with a clean cloth or paper towels.

Pay particular attention to the welding bar and its connectors (indicated by the arrow in the following figure) throughout the entire procedure described above for decontaminating and cleaning the cavity.

The welding bar can be removed to facilitate the cleaning of the cavity. Remember to reinsert it properly at the end of the procedure by pressing it downward into its lodging.

Check to make sure that the welding bar has been reinserted properly before using the unit again.

4) External surfaces

- Spray the external surfaces of the TissueSAFE plus unit with a disinfectant solution and let stand.
- Do not spray the unit's rear panel.
- Clean the surface with a water-soaked cloth to remove the disinfectant.
- Use an appropriate detergent to clean the surfaces.
- Clean the surface with a water-soaked cloth to remove the detergent.
- Dry the surface with a clean cloth or paper towels.

Do not spray directly into the printer, the barcode reader, the USB ports, or the slot for the temperature monitoring system (if present).

5) Final phases

- Reinsert the two spacers into the cavity.
- Wash your gloves with disinfectant, water, soap, and water again before removing them.
- Wash your hands thoroughly.
- Remove any remaining personal protective equipment.
- Make sure that the contaminated material is packaged and disposed of in an approved incinerator in accordance with the current local and national regulations.

At the end of the decontamination/cleaning operations, wait at least 5 minutes before using the unit or connecting it to its electrical power supply.

10. TRANSPORTATION AND STORAGE

The device, its components and accessories can be transported by rail, road, river and air transport (in covered vehicles).

- Milestone S.r.l. recommends to involve a specialized transport company to transport, install or possible handling of the device.
- Condensation water may form in the instrument if there is an extreme difference in temperature between the storage location and the installation site and if air humidity is high at the same time. If this is the case, wait at least two hours before switching on the instrument.

The fastening of a transport packaging in vehicles and the transportation of products must be carried out in accordance with the rules effective for the respective means of transport.

Arrangement and fastening of transport boxes should ensure their stable position and absence of displacement during transportation.

- The device may be transported only in an upright position.

Transport conditions: at -20°C...+60 °C. Relative humidity: up to 50% (at 40 °C), non-condensing.

Conditions of storage in transport packaging at the manufacturer's (user's) warehouses must comply with the following conditions: at -20°C...+60 °C. Relative humidity: up to 50% (at 40 °C), non-condensing. Rolls of printer paper should be stored separately at +10°C... +40 °C. Avoid direct exposure to sunlight.

It is not allowed to store products in a room where volatile liquids and substances that can cause corrosion keep.

After exposure to temperature and humidity during transportation and storage under the conditions specified in this technical file, the device remains fully functional.

11. MANUFACTURER WARRANTY

The manufacturer warrants that the quality of the unit meets the requirements of this Operator manual.

The device is assembled for a specific order, so there is no warranted storage period.

The warranty period for the device operation is 1 year from the date of sale, confirmed by the stamp and signature of the seller.

During the entire warranty period, the manufacturer undertakes to eliminate the defects free of charge or replace the defective device, if the damage is not related to the possible cases described below.

Possible cases that are not covered by this warranty:

- 1) Man-made damage or damage caused by misuse.
- 2) Damage caused by improper handling during transportation and storage of the device.
- 3) Damage caused by uncontrolled natural forces such as earthquake, fire or war.
- 4) The environment (operation conditions) in which the device is used does not meet the requirements specified in this Operator manual.

- 5) Damage caused by using the device with a power supply that does not meet the specifications in this Operator manual.

- 6) Damage caused by maintenance performed by personnel not authorized by the manufacturer Milestone S.r.l., Italy

- 7) Malfunction of the device, the registration number of which is not clearly legible.

Repair of the device is carried out only by Milestone S.r.l. It is strictly forbidden to open the device without an authorized representative of the manufacturer. The manufacturer warrants the continued safety of the equipment after maintenance and checks. The fuses integrated in the device are overcurrent fuses, and they are not subject to replacement by users. If you need to replace the fuse, please notify the Milestone service engineer.

If signs of unauthorized opening and mechanical damage are found, the manufacturer disclaims responsibility for the warranty obligations.

Service life of the device is 10 years.

Service life of battery in the temperature monitoring system is 2 years.

For maintenance and repair of the unit please contact the manufacturer or its authorized representative.

Manufacturer: Milestone S.r.l.

address: Via Fatebenefratelli, 1/5, IT-24010 Sorisole (BG), Italy

Web: www.milestonemed.com

tel.: +39 035 572735

fax: +39 035 575498

Authorized Representative in the Russian Federation: BioLine, LLC

Legal address: Russia, 197022, Saint Petersburg, inner-city Municipality Aptekarsky Ostrov, 23 Professora Popova Str., Let. E, premises 3H/5H, room 208

Postal address: 197022, Saint Petersburg, 23 Professora Popova Str., Let. E

tel.: +7(812) 320-49-49

e-mail: service@bioline.ru

Hot line: 8-800-555-49-40.


A. APPENDIX A – PROPER OPERATING METHODS

A.1. HOW TO POSITION A BAG IN THE TISSUESAFE PLUS UNIT

The bag containing the sample must be inserted with the opening facing the operator.

The plastic bag must be positioned so that the lines printed on the bag itself coincide with the welding bar inside the cavity, as shown in the following figure. Always position the uppermost line at the level of the welding bar.

A.2. HOW TO INSERT A SAMPLE INTO A VACUUM BAG

There are two possible ways to properly insert a sample into a bag:

1. If two operators are available: one can hold the bag open and the other can insert the sample into the bag without getting the bag's internal and external upper edges dirty.
2. If only one operator is available: turn down the upper external edges of the bag as shown in the following image before attempting to insert the sample. This will allow the operator to insert the sample without getting the bag's internal and external upper edges dirty.

A.3. HOW TO PROPERLY SEAL A VACUUM BAG

The instructions provided above must be followed in order to properly seal a vacuum bag, as the internal part of the bag must be clean and dry when the seal is applied. It is also extremely important to properly position the bag on the welding bar, as explained in this section.

At the end of a program, it is recommended to check whether the seal is uniform along the entire surface of the bag.

If the weld is not correct:

- Open the bag with scissors.
- Clean the upper inside part of the bag, as described above.
- Repeat the vacuum program.

B. APPENDIX B – VACUUM PROGRAMS GUIDELINES

Organs List:

P1	P2	P3	P4
Kidney (P3/P4 if with cyst)	Placenta	Lung	Cysts >10 cm in diameter
Breast	Prostate	Brain	
Uterus	Intestine	Lymph node	
Liver (P3/P4 if with cyst)	Testicle	Closed colon	
Muscle	Appendix	Gallbladder	
Thyroid	Bone (not sharp)**	Eye	
Thymus	Stomach / esophagus	Cysts <10cm in diameter	
Skin	Ovary	Fetus	
Tongue	Tonsil	Scraping/TURP	
Spleen			
Open colon*			
Pancreas			
Bladder			
Heart			

* if the colon is clamped, use program P3.

** in the case of bones with sharp edges, these must be covered with gauze in order to protect the vacuum bag.


If the nature of the specimen is uncertain, use the program P3.

This guide is intended for specimens directly packaged with Milestone vacuum sealed plastic bags.

C. APPENDIX C – TRANSPORTATION GUIDELINES

If the time necessary to transport the samples to Pathology is:

- **Less than 30 minutes:** place the sealed samples in the regigerator box with cool packs at +4°C.
- **Less than 1 hour:** place the samples in the refrigerator next to the TissueSAFE plus unit for at least 2 hours (pre-cooling procedure), and then move them into the regigerator box with cool packs at +4°C.
- **Greater than 3 hours:** use the mobile refrigerator unit connected to an external power supply (12/24Vdc, 230V~).

The controlled-temperature system designed by Milestone s.r.l. guarantees the perfect preservation of the sample's morphology; antigenicity, DNA/RNA, and proteins are preserved in an optimal manner up to 72 hours (depending on the initial conditions and the type of sample).


The total transport time must never exceed 72 hours in any of these three cases (including the time spent in the Operating Room or in the Histology Laboratory's refrigerator). In order to ensure the proper preservation of the samples' morphology, they must be stored and preserved at a temperature of +4°C±2°C.

D. APPENDIX D – PRINTER AND BARCODE READER

D.1. How to use the printer


A printer is included with the TissueSAFE plus unit in order to provide complete documentation regarding the program executed by printing an adhesive label that identifies the sample in the bag.

This printer allows the user to print 1D or 2D barcodes, based on the settings applied by the technical support personnel when configuring the software at the factory. By default, the instrument prints barcodes in 1D format, using the CODE 128B. The following figure shows the various formats that can be used by the printer, which can only be activated by the Technical Support Personnel.

A kit with 10 rolls for CODE 68331/A labels is included.

The printer consists of an ABS body with a front door, which provide access to the roll of paper and the printing mechanism. The console on the front (see the following figure) includes a FEED button (3) and an OPEN button for opening the paper roll compartment. The OPEN button is made from transparent plastic material, thus allowing it to also function as an LED STATUS indicator light (2).

- FEED button (3): when pressed, the printer feeds the paper through the mechanism.
- OPEN button (1): press this button to access the paper roll compartment.
- LED STATUS indicator (2): indicates the printer's hardware status. In the event of a malfunction, the flashing interval changes based on the indications provided in the following table.

STATUS LED	COLOR	DESCRIPTION	
	OFF	Printer OFF	
	GREEN	Printer ON: no error	
	GREEN	COMMUNICATION STATUS	
		Flashing	Description
		Slow	Paper end
		Fast	Heading over temperature or power supply voltage incorrect

The printer comes with a roll of paper already inserted, and another 9 rolls included separately within the packaging.

Once a roll is finished it must be replaced if you want to continue printing. The following instructions explain how to replace a roll of paper.

If you want to deactivate the printing function or if no roll of paper is available and therefore the printer cannot print for a certain period of time, the printing function can be deactivated in the "Owner Settings" (chapter 3.4).

D.1.1. Storing the rolls of paper


Temperature: +10°C...+40°C.

Avoid exposure to direct sunlight.

D.1.2. Changing the roll of paper

In order to change the roll of paper do the following:

- Open the printer's plastic cover by pressing the OPEN button.


- Lower the plastic cover until it is fully open.


Before inserting the new roll, cut off the first 20 cm of paper, as they might be sticky.

- Position the roll of paper inside the compartment, respecting the direction of rotation indicated.

- Pull the paper until it comes out of the compartment and close the door.
- Push on the plastic cover to lock it in place.
- Tear off the excess paper using the toothed edge.

D.1.3. Printer troubleshooting

It is necessary to use the toothed edge to tear off the paper.

If the paper does not advance properly, check to make sure that the roll is inserted correctly by following the operations described in the previous chapter.

It is also recommended to use the FEED button after each printing operation. This button causes the paper to advance slightly, thus preventing it from jamming the next time the printer is used.

 MILESTONE HELPING PATIENTS	 Milestone srl Via Fatebenefratelli, 1/5 - 24010 Sorisole (BG), Italy Tel. +39 035 4128264 - Fax +39 035 575498	marketing@milestonemedsrl.com www.milestonemed.com
---	--	--

CAMERA DI COMMERCIO, INDUSTRIA, ARTIGIANATO E AGRICOLTURA
CHAMBER OF COMMERCE, INDUSTRY, HANDICRAFT AND AGRICULTURE
BERGAMO

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We certify that the signer is a duly authorized signatory on behalf of the Company for any documents and acts used in trade with foreign countries.

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
Full list	Events Log					List length: 81	
Number	Type	Description	Date	Time	User	File	P
01	ACT	Start status Tef. samples 0	20-Oct-15	14:40:59	Admin (0001)	1	PS
02	ACT	Stop status Tef. samples 0	20-Oct-15	15:54:59	Admin (0001)	1	PS
03	ACT	Start status Tef. samples 0	20-Oct-15	16:28:59	Admin (0001)	1	PS
04	ACT	Stop status Tef. samples 0	20-Oct-15	17:28:59	Admin (0001)	1	PS
05	ACT	Start status Tef. samples 0	20-Oct-15	18:28:59	Admin (0001)	1	PS
06	ACT	Stop status Tef. samples 0	20-Oct-15	19:28:59	Admin (0001)	1	PS
07	ACT	Start status Tef. samples 0	20-Oct-15	20:28:59	Admin (0001)	1	PS
08	ACT	Stop status Tef. samples 0	20-Oct-15	21:28:59	Admin (0001)	1	PS
09	ACT	Start status Tef. samples 0	20-Oct-15	22:28:59	Admin (0001)	1	PS
10	ACT	Stop status Tef. samples 0	20-Oct-15	23:28:59	Admin (0001)	1	PS
11	ACT	Start status Tef. samples 0	20-Oct-15	00:28:59	Admin (0001)	1	PS
12	ACT	Stop status Tef. samples 0	20-Oct-15	01:28:59	Admin (0001)	1	PS
13	ACT	Start status Tef. samples 0	20-Oct-15	02:28:59	Admin (0001)	1	PS
14	ACT	Stop status Tef. samples 0	20-Oct-15	03:28:59	Admin (0001)	1	PS
15	ACT	Start status Tef. samples 0	20-Oct-15	04:28:59	Admin (0001)	1	PS
16	ACT	Stop status Tef. samples 0	20-Oct-15	05:28:59	Admin (0001)	1	PS
17	ACT	Start status Tef. samples 0	20-Oct-15	06:28:59	Admin (0001)	1	PS
18	ACT	Stop status Tef. samples 0	20-Oct-15	07:28:59	Admin (0001)	1	PS
19	ACT	Start status Tef. samples 0	20-Oct-15	08:28:59	Admin (0001)	1	PS
20	ACT	Stop status Tef. samples 0	20-Oct-15	09:28:59	Admin (0001)	1	PS
21	ACT	Start status Tef. samples 0	20-Oct-15	10:28:59	Admin (0001)	1	PS
22	ACT	Stop status Tef. samples 0	20-Oct-15	11:28:59	Admin (0001)	1	PS
23	ACT	Start status Tef. samples 0	20-Oct-15	12:28:59	Admin (0001)	1	PS
24	ACT	Stop status Tef. samples 0	20-Oct-15	13:28:59	Admin (0001)	1	PS
25	ACT	Start status Tef. samples 0	20-Oct-15	14:28:59	Admin (0001)	1	PS
26	ACT	Stop status Tef. samples 0	20-Oct-15	15:28:59	Admin (0001)	1	PS
27	ACT	Start status Tef. samples 0	20-Oct-15	16:28:59	Admin (0001)	1	PS
28	ACT	Stop status Tef. samples 0	20-Oct-15	17:28:59	Admin (0001)	1	PS
29	ACT	Start status Tef. samples 0	20-Oct-15	18:28:59	Admin (0001)	1	PS
30	ACT	Stop status Tef. samples 0	20-Oct-15	19:28:59	Admin (0001)	1	PS
31	ACT	Start status Tef. samples 0	20-Oct-15	20:28:59	Admin (0001)	1	PS
32	ACT	Stop status Tef. samples 0	20-Oct-15	21:28:59	Admin (0001)	1	PS
33	ACT	Start status Tef. samples 0	20-Oct-15	22:28:59	Admin (0001)	1	PS
34	ACT	Stop status Tef. samples 0	20-Oct-15	23:28:59	Admin (0001)	1	PS
35	ACT	Start status Tef. samples 0	20-Oct-15	00:28:59	Admin (0001)	1	PS
36	ACT	Stop status Tef. samples 0	20-Oct-15	01:28:59	Admin (0001)	1	PS
37	ACT	Start status Tef. samples 0	20-Oct-15	02:28:59	Admin (0001)	1	PS
38	ACT	Stop status Tef. samples 0	20-Oct-15	03:28:59	Admin (0001)	1	PS
39	ACT	Start status Tef. samples 0	20-Oct-15	04:28:59	Admin (0001)	1	PS
40	ACT	Stop status Tef. samples 0	20-Oct-15	05:28:59	Admin (0001)	1	PS
41	ACT	Start status Tef. samples 0	20-Oct-15	06:28:59	Admin (0001)	1	PS
42	ACT	Stop status Tef. samples 0	20-Oct-15	07:28:59	Admin (0001)	1	PS
43	ACT	Start status Tef. samples 0	20-Oct-15	08:28:59	Admin (0001)	1	PS
44	ACT	Stop status Tef. samples 0	20-Oct-15	09:28:59	Admin (0001)	1	PS
45	ACT	Start status Tef. samples 0	20-Oct-15	10:28:59	Admin (0001)	1	PS
46	ACT	Stop status Tef. samples 0	20-Oct-15	11:28:59	Admin (0001)	1	PS
47	ACT	Start status Tef. samples 0	20-Oct-15	12:28:59	Admin (0001)	1	PS
48	ACT	Stop status Tef. samples 0	20-Oct-15	13:28:59	Admin (0001)	1	PS
49	ACT	Start status Tef. samples 0	20-Oct-15	14:28:59	Admin (0001)	1	PS
50	ACT	Stop status Tef. samples 0	20-Oct-15	15:28:59	Admin (0001)	1	PS
51	ACT	Start status Tef. samples 0	20-Oct-15	16:28:59	Admin (0001)	1	PS
52	ACT	Stop status Tef. samples 0	20-Oct-15	17:28:59	Admin (0001)	1	PS
53	ACT	Start status Tef. samples 0	20-Oct-15	18:28:59	Admin (0001)	1	PS
54	ACT	Stop status Tef. samples 0	20-Oct-15	19:28:59	Admin (0001)	1	PS
55	ACT	Start status Tef. samples 0	20-Oct-15	20:28:59	Admin (0001)	1	PS
56	ACT	Stop status Tef. samples 0	20-Oct-15	21:28:59	Admin (0001)	1	PS
57	ACT	Start status Tef. samples 0	20-Oct-15	22:28:59	Admin (0001)	1	PS
58	ACT	Stop status Tef. samples 0	20-Oct-15	23:28:59	Admin (0001)	1	PS
59	ACT	Start status Tef. samples 0	20-Oct-15	00:28:59	Admin (0001)	1	PS
60	ACT	Stop status Tef. samples 0	20-Oct-15	01:28:59	Admin (0001)	1	PS
61	ACT	Start status Tef. samples 0	20-Oct-15	02:28:59	Admin (0001)	1	PS
62	ACT	Stop status Tef. samples 0	20-Oct-15	03:28:59	Admin (0001)	1	PS
63	ACT	Start status Tef. samples 0	20-Oct-15	04:28:59	Admin (0001)	1	PS
64	ACT	Stop status Tef. samples 0	20-Oct-15	05:28:59	Admin (0001)	1	PS
65	ACT	Start status Tef. samples 0	20-Oct-15	06:28:59	Admin (0001)	1	PS
66	ACT	Stop status Tef. samples 0	20-Oct-15	07:28:59	Admin (0001)	1	PS
67	ACT	Start status Tef. samples 0	20-Oct-15	08:28:59	Admin (0001)	1	PS
68	ACT	Stop status Tef. samples 0	20-Oct-15	09:28:59	Admin (0001)	1	PS
69	ACT	Start status Tef. samples 0	20-Oct-15	10:28:59	Admin (0001)	1	PS
70	ACT	Stop status Tef. samples 0	20-Oct-15	11:28:59	Admin (0001)	1	PS
71	ACT	Start status Tef. samples 0	20-Oct-15	12:28:59	Admin (0001)	1	PS
72	ACT	Stop status Tef. samples 0	20-Oct-15	13:28:59	Admin (0001)	1	PS
73	ACT	Start status Tef. samples 0	20-Oct-15	14:28:59	Admin (0001)	1	PS
74	ACT	Stop status Tef. samples 0	20-Oct-15	15:28:59	Admin (0001)	1	PS
75	ACT	Start status Tef. samples 0	20-Oct-15	16:28:59	Admin (0001)	1	PS
76	ACT	Stop status Tef. samples 0	20-Oct-15	17:28:59	Admin (0001)	1	PS
77	ACT	Start status Tef. samples 0	20-Oct-15	18:28:59	Admin (0001)	1	PS
78	ACT	Stop status Tef. samples 0	20-Oct-15	19:28:59	Admin (0001)	1	PS
79	ACT	Start status Tef. samples 0	20-Oct-15	20:28:59	Admin (0001)	1	PS
80	ACT	Stop status Tef. samples 0	20-Oct-15	21:28:59	Admin (0001)	1	PS
81	ACT	Start status Tef. samples 0	20-Oct-15	22:28:59	Admin (0001)	1	PS

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Unit name: SediSAFE	
Unit name	TissueSAFE plus
Software Version	15.11
Serial number	125199002
Processing name	P1
Category	VACUUM ONLY
User	Admin
Start date	2015/12/09, at 11:01
Specimen weight	30 g
Outcome	completed with success
Bag ID	387912523028372A
Bag size	SelfStanding LARGE

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Saved
 LogEvents.tissuesafepius.s3db

➔

SUBJECT	TYPE	DESCRIPTION	DATE	TIME	STATUS	CTD	LOFF	LOON	SW
1	Aut	Registration	11 Nov 15	09:22:58	---	1	1	0	1058
2	Aut	User started	11 Nov 15	09:23:11	Service (1)	1	0	0	1060
3	Aut	Access to T	11 Nov 15	09:23:22	Service (1)	1	0	0	1062
4	Aut	Access to T	11 Nov 15	09:24:26	Service (1)	1	0	0	110
5	Aut	Set back to 1	11 Nov 15	09:24:48	Service (1)	1	0	0	112
6	Aut	Registration	11 Nov 15	09:24:51	Service (1)	1	0	0	114
7	Aut	Registration	11 Nov 15	09:24:54	Service (1)	1	0	0	116
8	Aut	Registration	11 Nov 15	09:24:56	Service (1)	1	0	0	118
9	Aut	Registration	11 Nov 15	09:24:58	Service (1)	1	0	0	120
10	Aut	Registration	11 Nov 15	09:25:00	Service (1)	1	0	0	122
11	Aut	Registration	11 Nov 15	09:25:02	Service (1)	1	0	0	124
12	Aut	Registration	11 Nov 15	09:25:04	Service (1)	1	0	0	126
13	Aut	Registration	11 Nov 15	09:25:06	Service (1)	1	0	0	128
14	Aut	Registration	11 Nov 15	09:25:08	Service (1)	1	0	0	130
15	Aut	Registration	11 Nov 15	09:25:10	Service (1)	1	0	0	132
16	Aut	Registration	11 Nov 15	09:25:12	Service (1)	1	0	0	134
17	Aut	Registration	11 Nov 15	09:25:14	Service (1)	1	0	0	136
18	Aut	Registration	11 Nov 15	09:25:16	Service (1)	1	0	0	138
19	Aut	Registration	11 Nov 15	09:25:18	Service (1)	1	0	0	140
20	Aut	Registration	11 Nov 15	09:25:20	Service (1)	1	0	0	142
21	Aut	Registration	11 Nov 15	09:25:22	Service (1)	1	0	0	144
22	Aut	Registration	11 Nov 15	09:25:24	Service (1)	1	0	0	146
23	Aut	Registration	11 Nov 15	09:25:26	Service (1)	1	0	0	148
24	Aut	Registration	11 Nov 15	09:25:28	Service (1)	1	0	0	150
25	Aut	Registration	11 Nov 15	09:25:30	Service (1)	1	0	0	152
26	Aut	Registration	11 Nov 15	09:25:32	Service (1)	1	0	0	154
27	Aut	Registration	11 Nov 15	09:25:34	Service (1)	1	0	0	156
28	Aut	Registration	11 Nov 15	09:25:36	Service (1)	1	0	0	158
29	Aut	Registration	11 Nov 15	09:25:38	Service (1)	1	0	0	160
30	Aut	Registration	11 Nov 15	09:25:40	Service (1)	1	0	0	162
31	Aut	Registration	11 Nov 15	09:25:42	Service (1)	1	0	0	164
32	Aut	Registration	11 Nov 15	09:25:44	Service (1)	1	0	0	166
33	Aut	Registration	11 Nov 15	09:25:46	Service (1)	1	0	0	168
34	Aut	Registration	11 Nov 15	09:25:48	Service (1)	1	0	0	170
35	Aut	Registration	11 Nov 15	09:25:50	Service (1)	1	0	0	172
36	Aut	Registration	11 Nov 15	09:25:52	Service (1)	1	0	0	174
37	Aut	Registration	11 Nov 15	09:25:54	Service (1)	1	0	0	176
38	Aut	Registration	11 Nov 15	09:25:56	Service (1)	1	0	0	178
39	Aut	Registration	11 Nov 15	09:25:58	Service (1)	1	0	0	180
40	Aut	Registration	11 Nov 15	09:26:00	Service (1)	1	0	0	182
41	Aut	Registration	11 Nov 15	09:26:02	Service (1)	1	0	0	184
42	Aut	Registration	11 Nov 15	09:26:04	Service (1)	1	0	0	186
43	Aut	Registration	11 Nov 15	09:26:06	Service (1)	1	0	0	188
44	Aut	Registration	11 Nov 15	09:26:08	Service (1)	1	0	0	190
45	Aut	Registration	11 Nov 15	09:26:10	Service (1)	1	0	0	192
46	Aut	Registration	11 Nov 15	09:26:12	Service (1)	1	0	0	194
47	Aut	Registration	11 Nov 15	09:26:14	Service (1)	1	0	0	196
48	Aut	Registration	11 Nov 15	09:26:16	Service (1)	1	0	0	198
49	Aut	Registration	11 Nov 15	09:26:18	Service (1)	1	0	0	200
50	Aut	Registration	11 Nov 15	09:26:20	Service (1)	1	0	0	202
51	Aut	Registration	11 Nov 15	09:26:22	Service (1)	1	0	0	204
52	Aut	Registration	11 Nov 15	09:26:24	Service (1)	1	0	0	206
53	Aut	Registration	11 Nov 15	09:26:26	Service (1)	1	0	0	208
54	Aut	Registration	11 Nov 15	09:26:28	Service (1)	1	0	0	210
55	Aut	Registration	11 Nov 15	09:26:30	Service (1)	1	0	0	212
56	Aut	Registration	11 Nov 15	09:26:32	Service (1)	1	0	0	214
57	Aut	Registration	11 Nov 15	09:26:34	Service (1)	1	0	0	216
58	Aut	Registration	11 Nov 15	09:26:36	Service (1)	1	0	0	218
59	Aut	Registration	11 Nov 15	09:26:38	Service (1)	1	0	0	220
60	Aut	Registration	11 Nov 15	09:26:40	Service (1)	1	0	0	222
61	Aut	Registration	11 Nov 15	09:26:42	Service (1)	1	0	0	224
62	Aut	Registration	11 Nov 15	09:26:44	Service (1)	1	0	0	226
63	Aut	Registration	11 Nov 15	09:26:46	Service (1)	1	0	0	228
64	Aut	Registration	11 Nov 15	09:26:48	Service (1)	1	0	0	230
65	Aut	Registration	11 Nov 15	09:26:50	Service (1)	1	0	0	232
66	Aut	Registration	11 Nov 15	09:26:52	Service (1)	1	0	0	234
67	Aut	Registration	11 Nov 15	09:26:54	Service (1)	1	0	0	236
68	Aut	Registration	11 Nov 15	09:26:56	Service (1)	1	0	0	238
69	Aut	Registration	11 Nov 15	09:26:58	Service (1)	1	0	0	240
70	Aut	Registration	11 Nov 15	09:27:00	Service (1)	1	0	0	242
71	Aut	Registration	11 Nov 15	09:27:02	Service (1)	1	0	0	244
72	Aut	Registration	11 Nov 15	09:27:04	Service (1)	1	0	0	246
73	Aut	Registration	11 Nov 15	09:27:06	Service (1)	1	0	0	248
74	Aut	Registration	11 Nov 15	09:27:08	Service (1)	1	0	0	250
75	Aut	Registration	11 Nov 15	09:27:10	Service (1)	1	0	0	252
76	Aut	Registration	11 Nov 15	09:27:12	Service (1)	1	0	0	254
77	Aut	Registration	11 Nov 15	09:27:14	Service (1)	1	0	0	256
78	Aut	Registration	11 Nov 15	09:27:16	Service (1)	1	0	0	258
79	Aut	Registration	11 Nov 15	09:27:18	Service (1)	1	0	0	260
80	Aut	Registration	11 Nov 15	09:27:20	Service (1)	1	0	0	262
81	Aut	Registration	11 Nov 15	09:27:22	Service (1)	1	0	0	264
82	Aut	Registration	11 Nov 15	09:27:24	Service (1)	1	0	0	266
83	Aut	Registration	11 Nov 15	09:27:26	Service (1)	1	0	0	268
84	Aut	Registration	11 Nov 15	09:27:28	Service (1)	1	0	0	270
85	Aut	Registration	11 Nov 15	09:27:30	Service (1)	1	0	0	272
86	Aut	Registration	11 Nov 15	09:27:32	Service (1)	1	0	0	274
87	Aut	Registration	11 Nov 15	09:27:34	Service (1)	1	0	0	276
88	Aut	Registration	11 Nov 15	09:27:36	Service (1)	1	0	0	278
89	Aut	Registration	11 Nov 15	09:27:38	Service (1)	1	0	0	280
90	Aut	Registration	11 Nov 15	09:27:40	Service (1)	1	0	0	282
91	Aut	Registration	11 Nov 15	09:27:42	Service (1)	1	0	0	284
92	Aut	Registration	11 Nov 15	09:27:44	Service (1)	1	0	0	286
93	Aut	Registration	11 Nov 15	09:27:46	Service (1)	1	0	0	288
94	Aut	Registration	11 Nov 15	09:27:48	Service (1)	1	0	0	290
95	Aut	Registration	11 Nov 15	09:27:50	Service (1)	1	0	0	292
96	Aut	Registration	11 Nov 15	09:27:52	Service (1)	1	0	0	294
97	Aut	Registration	11 Nov 15	09:27:54	Service (1)	1	0	0	296
98	Aut	Registration	11 Nov 15	09:27:56	Service (1)	1	0	0	298
99	Aut	Registration	11 Nov 15	09:27:58	Service (1)	1	0	0	300
100	Aut	Registration	11 Nov 15	09:28:00	Service (1)	1	0	0	302

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