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## CERTIFICATE

中国国际贸易促进委员会暨中国国际商会  
China Council for the Promotion of International Trade is China Chamber of International Commerce

01857376

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China Council for the Promotion of International Trade  
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HUQIU IMAGING (SUZHOU) CO., LTD on the annexed DOCUMENT is genuine. The Russian translation of the attached document is consistent with the English version.

China Council for the Promotion  
of International Trade

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Authorized

Signature:

Sun Danfeng

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(Date: Oct. 13, 2023)

查询网址: Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

**Approved**

Huqiu Imaging (Suzhou) Co., Ltd  
Manufacturer

No.19 Tianmushan Road, SND, Suzhou 215163, Jiangsu, China  
Legal Address

Sales Director  
Occupation

Allen Wang  
Surname, Name

29.09.2023  
Date DD.MM.YYYY

Signature

**INSTRUCTIONS FOR USE  
OF THE MEDICAL DEVICE**  
DEVICE FOR PRINTING DIGITAL DIAGNOSTIC MEDICAL IMAGES OF MODELS:  
**MS-460DY, MS-760DY**

**2023**

Document version v2.20.09.23

## CONTENTS

| NAME                                                                                                                    | PAGE |
|-------------------------------------------------------------------------------------------------------------------------|------|
| ABOUT SAFETY NOTES IN THIS DOCUMENT                                                                                     | 5    |
| DEVICE FOR PRINTING DIGITAL DIAGNOSTIC MEDICAL IMAGES OF MODELS:<br>MS-460DY, MS-760DY.                                 | 6    |
| INFORMATION ABOUT THE MANUFACTURER OF THE MEDICAL DEVICE (OR) ITS<br>AUTHORIZED REPRESENTATIVE                          | 7    |
| APPLICATION AREA                                                                                                        | 8    |
| PURPOSE                                                                                                                 | 9    |
| INDICATIONS AND CONTRAINDICATIONS FOR THE USE OF MEDICAL DEVICE                                                         | 10   |
| POSSIBLE SIDE EFFECTS                                                                                                   | 11   |
| DESCRIPTION OF THE PRINCIPLES ON WHICH THE OPERATION OF THE MEDICAL<br>DEVICE IS BASED                                  | 12   |
| FILM TYPES                                                                                                              | 13   |
| SPECIFICATIONS / SCOPE OF SUPPLY                                                                                        | 14   |
| PRODUCT CLAIMS                                                                                                          | 17   |
| INTRODUCTION                                                                                                            | 18   |
| INTENDED CONSUMERS AND USERS                                                                                            | 19   |
| FUNCTIONS AND FEATURES                                                                                                  | 20   |
| EQUIPMENT CLASSIFICATION                                                                                                | 21   |
| OPERATING RISKS AND PRECAUTIONS                                                                                         | 22   |
| SYMBOLS AND CONVENTIONS                                                                                                 | 23   |
| INSTALLATION                                                                                                            | 25   |
| SPACE REQUIREMENTS FOR THE ROOM                                                                                         | 26   |
| TRAINING                                                                                                                | 27   |
| LIST OF APPLICABLE STANDARDS                                                                                            | 28   |
| INTERACTION WITH EXTERNAL SYSTEMS                                                                                       | 29   |
| MARKING                                                                                                                 | 30   |
| INSTRUCTIONS FOR USING POWER CABLES                                                                                     | 31   |
| PATIENT DATA SECURITY                                                                                                   | 32   |
| OPERATING MODES                                                                                                         | 34   |
| LOCAL USER INTERFACE                                                                                                    | 35   |
| TURN ON THE DEVICE FOR PRINTING                                                                                         | 39   |
| COOLING DEVICE FOR PRINTING                                                                                             | 41   |
| TURN OFF THE PRINT DEVICE                                                                                               | 42   |
| BASIC OPERATIONS (USER MODE)                                                                                            | 43   |
| PRINT TASK MANAGEMENT                                                                                                   | 44   |
| PAUSE PRINT TASK                                                                                                        | 45   |
| VIEW THE NUMBER OF FILMS                                                                                                | 46   |
| DELETE PRINTING TASK                                                                                                    | 48   |
| PRINT SIZE                                                                                                              | 49   |
| FILM LOADING                                                                                                            | 50   |
| PROPER POSITIONING OF THE FILM IN THE LOADING TRAY                                                                      | 53   |
| MENU STRUCTURE                                                                                                          | 54   |
| QUALITY CONTROL                                                                                                         | 55   |
| TEST IMAGE FOR QUALITY CONTROL                                                                                          | 56   |
| CHECKING THE ALLOWABLE SPATIAL RESOLUTION, THE LEVEL OF INTERFERENCE<br>AND THE VISIBILITY OF DETAILS WITH LOW CONTRAST | 57   |
| DETERMINATION OF THE BASIC VALUES OF THE GEOMETRIC PARAMETERS OF THE<br>IMAGE                                           | 58   |
| SCHEDULE OF PREVENTIVE CARE ACTIVITIES                                                                                  | 59   |
| SAFETY INSTRUCTIONS                                                                                                     | 60   |
| CLEANING AND DISINFECTION                                                                                               | 61   |
| REQUIREMENTS FOR THE TECHNICAL SERVICE AND REPAIR OF A MEDICAL DEVICE                                                   | 62   |
| CHECKING AND CONTROL OF THE SYSTEM BY PERSONNEL USING THE DEVICE                                                        | 63   |

Document version v2.20.09.23

|                                                                                  |    |
|----------------------------------------------------------------------------------|----|
| TECHNICAL SERVICE OF THE SYSTEM AT THE CONSUMER                                  | 64 |
| MOUNTING AND INSTALLATION REQUIREMENT                                            | 65 |
| THE PROCEDURE AND CONDITIONS FOR THE RECYCLING OR DISPOSAL OF THE MEDICAL DEVICE | 66 |
| REQUIREMENTS FOR ENVIRONMENTAL PROTECTION WHEN USING THE MEDICAL DEVICE          | 67 |
| WARRANTY OBLIGATIONS AND RECLAMATION OF THE MEDICAL DEVICE                       | 68 |
| HIGH-FREQUENCY RADIATION AND PROTECTION INFORMATION                              | 69 |
| RF IMMUNITY OF WIRELESS COMMUNICATION EQUIPMENT                                  | 74 |
| EMC (ELECTROMAGNETIC CAPABILITY) PRECAUTIONS                                     | 75 |
| CABLES AND ACCESSORIES                                                           | 77 |
| APPENDIX A (CODES AND TYPES OF ERRORS)                                           | 78 |

The present INSTRUCTIONS FOR USE applies to Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY and contains information about the device, safety precautions for operation, basic operations, and care and maintenance.

Read the INSTRUCTIONS FOR USE carefully before you start working with Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY. Pay special attention to warnings and cautions!

INSTRUCTIONS FOR USE may be supplied on electronic media. On request, a paper version is provided.

## OUT SAFETY NOTES IN THIS DOCUMENT

The following are examples of representing prescriptions of the type «Warning», «Attention», «Instruction» and «Note» on the pages of this document. The text of the examples explains the meaning of the respective warning / prescriptive block.

|                                                                                     |                                                                                                                                                                                                                                                                         |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <b>DANGER:</b><br>A «Danger» statement indicates a situation of direct, imminent danger of serious injury to the operator, engineer, patient, or others.                                                                                                                |
|    | <b>WARNING:</b><br>A «Warning» statement indicates a situation that could result in serious injury to the operator, engineer, patient, or others.                                                                                                                       |
|    | <b>ATTENTION:</b><br>An «Attention» statement indicates a situation that could result in minor injury to the operator, engineer, patient, or others.                                                                                                                    |
|    | <b>INSTRUCTION:</b><br>«Instruction» statement contain instructions that, if not followed, could result in damage to the equipment mentioned in the INSTRUCTIONS FOR USE or other equipment or property, and may result in environmental pollution.                     |
|   | <b>PROHIBITION (DO NOT):</b><br>A «Prohibition (Do not)» statement contains instructions that, if not followed, could cause damage to the equipment referenced in the INSTRUCTIONS FOR USE or other equipment or property, and could result in environmental pollution. |
|  | <b>NOTE:</b><br>«Note» contain recommendations or clarification of special points. The «Note» contains no instructions.                                                                                                                                                 |

**VICE FOR PRINTING DIGITAL DIAGNOSTIC MEDICAL IMAGES OF MODELS:  
MS-460DY, MS-760DY**

---

- Information about the manufacturer of the MEDICAL DEVICE and (or) its authorized representative
- Application area
- Purpose
- Film types
- Specifications
- Product claims

Document version v2.20.09.23

**INFORMATION ABOUT THE MANUFACTURER OF THE MEDICAL DEVICE (OR) ITS  
AUTHORIZED REPRESENTATIVE**

---

**Manufacturer's name:** Huqiu Imaging (Suzhou)Co.,Ltd.

**Manufacturer's address:** No.19 Tianmushan Road, SND, Suzhou 215163, Jiangsu, China

**Developer:** Huqiu Imaging (Suzhou)Co., Ltd

**Place of production:** No.19 Tianmushan Road, SND, Suzhou 215163, Jiangsu, China

**Authorized representative of the manufacturer in the territory of the Russian Federation:**

Medcorp-MT, Stromynka str., building 19, housing. 2, Э 1 ПИВБ К32 PM3 107076 Moscow, Russia  
Tel. +7 (495) 969-61-47, medcorp@medcorp.ru

## PLICATION AREA

Device for printing digital diagnostic medical images of models (hereinafter referred to as Device for printing), of models: MS-460DY, MS-760DY can only be used for studies in general radiography. The MS-760DY can also be used for mammography studies.

Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY are designed for high productivity and can be used as central network devices for printing for printing medical images transmitted over the hospital network in DICOM format.

Image sources are images obtained from Computer Tomography Imaging, Magnetic Resonance Tomography, digital and computed radiography systems, fluoroscopy, PACS, nuclear medicine, mobile radiography, dentistry and other medical imaging systems.

**Table 1: Differences between the two Device for printing models**

| Device for printing digital diagnostic medical images,<br>model: MS-460DY | Device for printing digital diagnostic medical images,<br>model: MS-760DY |
|---------------------------------------------------------------------------|---------------------------------------------------------------------------|
| Maximum print resolution<br>320dpi (12px/mm)                              | Maximum print resolution<br>508dpi (20px/mm)                              |

## POSE

Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY is desktop digital thermal device designed for thermal printing of digital diagnostic medical images.

Image sources are images obtained from Computer Tomography Imaging, Magnetic Resonance Tomography, digital and computed radiography systems, fluoroscopy, PACS, nuclear medicine, mobile radiography, dentistry and other medical imaging systems.

## **INDICATIONS AND CONTRAINDICATIONS FOR THE USE OF THE MEDICAL DEVICE**

### **Indications**

Device for printing is intended for use in combination with the original Thermographic film for medical imaging systems «Dry Film MS-KX410 A/F», manufactured by Huqiu Imaging (Suzhou)Co., Ltd, for diagnosing normal or pathological human conditions in the field of radiography.

### **Contraindications**

When used for its intended purpose, proper operation and training of personnel, there are no side effects on the body and human health.

## **POSSIBLE SIDE EFFECTS**

---

No side effects have been identified.

It is recommended to strictly follow the manufacturer's instructions regarding the use of this MEDICAL DEVICE.

Document version v2.20.09.23

11

## **DESCRIPTION OF THE PRINCIPLES ON WHICH THE OPERATION OF THE MEDICAL DEVICE IS BASED**

---

In Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY direct thermal printing technology is implemented, which eliminates complex optical components and allows you to get black and white images with high diagnostic quality. The technology is based on receiving a signal from a source for obtaining or storing a radiographic image and converting this signal into an indelible X-ray image, which is printed on thermographic film for medical imaging systems with a resolution of 320 dpi/508 dpi. The MEDICAL DEVICE consists of an outer casing, a stand, a film entry mechanism, a film ejection mechanism, a thermal head, an embossing roller, a control system, a film cassette, a power supply, a network interface, and a USB interface. This MEDICAL DEVICE can send print tasks with diagnostic images obtained from various medical equipment used in CR, DR, Computer Tomography Imaging, MR, Magnetic Resonance Tomography that support the DICOM3.0 standard through the RJ45 network port.

## M TYPES

| Film name                                                              | Film base color  | Application         | Film format                                                                                                     | Maximum optical density |
|------------------------------------------------------------------------|------------------|---------------------|-----------------------------------------------------------------------------------------------------------------|-------------------------|
| Thermographic film for medical imaging systems «Dry Film MS-KX410 A/F» | Blue-transparent | General radiography | 20 x 25 cm (8x10 inches)<br>25 x 30 cm (10x12 inches)<br>28 x 35 cm (11x14 inches)<br>35 x 43 cm (14x17 inches) | $\geq 2,15$             |


The Device for printing has 2 input trays. Both input trays can use the following types and sizes of film:  
Top tray for 20 x 25 cm (8x10 inches), 25 x 30 cm (10x12 inches), 28 x 35 cm (11x14 inches)  
Bottom tray for 28 x 35 cm (11x14 inches), 35 x 43 cm (14x17 inches).

When loading a new pack of film, it reads the film IDs and automatically makes the necessary settings.

### Tray markings

*Note:* To change the film size, you must change the tray configuration

Appropriate marks are placed on the trays by service personnel when the device for printing is installed and indicate the type of new film that should be used when loading the tray.

## SPECIFICATIONS/ SCOPE OF SUPPLY

| Product description                     |                                                                                                                                                                                                                                     |                                                                                                                                                                             |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Model                                   | MS-460DY                                                                                                                                                                                                                            | MS-760DY                                                                                                                                                                    |
| Manufacturer                            | Huqiu Imaging (Suzhou)Co., Ltd                                                                                                                                                                                                      |                                                                                                                                                                             |
| Place of production                     | No.19 Tianmushan Road, SND, Suzhou 215163, Jiangsu, China                                                                                                                                                                           |                                                                                                                                                                             |
| Performance                             | 14x17 inches (35 x 43 cm) ≥ 70 sheets/h;<br>11x14 inches (28 x 35 cm) ≥ 75 sheets/h;<br>10x12 inches (25 x 30 cm) ≥ 80 sheets/h;<br>8x10 inches (20 x 25 cm) ≥ 90 sheets/h;                                                         | 14x17 inches (35 x 43 cm) ≥ 60 sheets/h;<br>11x14 inches (28 x 35 cm) ≥ 65 sheets/h;<br>10x12 inches (25 x 30 cm) ≥ 70 sheets/h;<br>8x10 inches (20 x 25 cm) ≥ 80 sheets/h; |
| Main memory/ROM (Read-only memory)      | At least 2 GB / 120 GB                                                                                                                                                                                                              |                                                                                                                                                                             |
| External storage device                 | USB flash drive                                                                                                                                                                                                                     |                                                                                                                                                                             |
| External dimensions                     | 740±5mm<br>566±5mm<br>445±5mm                                                                                                                                                                                                       |                                                                                                                                                                             |
| Weight                                  | 50 kg±5%                                                                                                                                                                                                                            | 50 kg±5%                                                                                                                                                                    |
| Number of film trays                    | 2                                                                                                                                                                                                                                   |                                                                                                                                                                             |
| Print Method                            | Direct Thermal                                                                                                                                                                                                                      |                                                                                                                                                                             |
| Film type                               | Thermographic film for medical imaging systems «Dry Film MS-KX410 A/F»                                                                                                                                                              |                                                                                                                                                                             |
| Version                                 | 8x10 inches (20x25 cm) (package contains 100 sheets);<br>10x12 inches (25x30 cm) (package contains 100 sheets);<br>11x14 inches (28x35 cm) (package contains 100 sheets);<br>14x17 inches (35x43 cm) (package contains 100 sheets). |                                                                                                                                                                             |
| Network protocols (TCP/IP services)     | HTTP                                                                                                                                                                                                                                |                                                                                                                                                                             |
| Interface                               | 10/100/1000 Base-T Ethernet (Twisted pair RJ-45)                                                                                                                                                                                    |                                                                                                                                                                             |
| Network protocol                        | Standard DICOM 3.0                                                                                                                                                                                                                  |                                                                                                                                                                             |
| Maximum print resolution                | 320dpi (12px/mm)                                                                                                                                                                                                                    | 508dpi (20px/mm)                                                                                                                                                            |
| Number of shades of gray for each pixel | 14 bit                                                                                                                                                                                                                              |                                                                                                                                                                             |
| Image quality                           | Automatic calibration with built-in densitometer                                                                                                                                                                                    |                                                                                                                                                                             |
| Protection level                        | Electric Shock - Class 1 (grounded)                                                                                                                                                                                                 |                                                                                                                                                                             |

|                                                                          |                                                                                                 |
|--------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------|
|                                                                          | Water penetration-IPX0                                                                          |
| Noise level                                                              | During operation - Max. 64 dBA<br>Standby - Max. 54 dBA                                         |
| Power supply system                                                      |                                                                                                 |
| Input voltage                                                            | 100-240VAC                                                                                      |
| Power frequency                                                          | 50-60 Hz                                                                                        |
| Maximum power consumption                                                | 600 VA±10%                                                                                      |
| Terms of Use<br>Temperature<br>Relative humidity<br>Atmospheric pressure | +15°C to +35°C<br>From 20% to 75% Note: Moistening of films is not allowed!<br>70 kPa - 106 kPa |
| Terms of Transportation                                                  | -10°C to +50°C<br>10% to 80%<br>50 kPa - 106 kPa                                                |
| Storage conditions                                                       | -10°C to +50°C<br>10% to 80%<br>50 kPa - 106 kPa                                                |
| Network cable                                                            | Plug IEC C13                                                                                    |
| Cable for connecting to a computer network                               | UTP; AWG 24                                                                                     |

| Image matrix - area of the diagnostic image MS-460DY |                      |        |                      |        |
|------------------------------------------------------|----------------------|--------|----------------------|--------|
| 20 x 25cm.(8x10")                                    | Dimension 20cm.(8")  |        | Dimension 25cm.(10") |        |
|                                                      | pixels               | mm     | pixels               | mm     |
|                                                      | 2376                 | 188,65 | 3070                 | 243,76 |
| 25 x 30cm.(10x12")                                   | Dimension 25cm.(10") |        | Dimension 30cm.(12") |        |
|                                                      | pixels               | mm     | pixels               | mm     |
|                                                      | 3070                 | 243,76 | 3653                 | 290,05 |
| 28 x 35cm.(11x14")                                   | Dimension 28cm.(11") |        | Dimension 35cm.(14") |        |
|                                                      | pixels               | mm     | pixels               | mm     |
|                                                      | 3348                 | 265,83 | 4358                 | 346,03 |
| 35 x 43cm.(14x17")                                   | Dimension 35cm.(14") |        | Dimension 43cm.(17") |        |
|                                                      | pixels               | mm     | pixels               | mm     |
|                                                      | 4358                 | 346,03 | 5232                 | 415,42 |

| Image matrix - area of the diagnostic image MS-760DY |                      |       |                      |     |
|------------------------------------------------------|----------------------|-------|----------------------|-----|
| 20 x 25cm.(8x10")                                    | Dimension 20cm.(8")  |       | Dimension 25cm.(10") |     |
|                                                      | pixels               | mm    | pixels               | mm  |
|                                                      | 3852                 | 192,6 | 4880                 | 244 |
| 25 x 30cm.(10x12")                                   | Dimension 25cm.(10") |       | Dimension 30cm.(12") |     |
|                                                      | pixels               | mm    | pixels               | mm  |
|                                                      | 4880                 | 244   | 5860                 | 293 |
| 28 x 35cm.(11x14")                                   | Dimension 28cm.(11") |       | Dimension 35cm.(14") |     |
|                                                      | pixels               | mm    | pixels               | mm  |

Document version v2.20.09.23

|                       |                      |       |                      |       |
|-----------------------|----------------------|-------|----------------------|-------|
|                       | 5376                 | 268,8 | 6922                 | 346,1 |
| 35 x<br>43cm.(14x17") | Dimension 35cm.(14") |       | Dimension 43cm.(17") |       |
|                       | pixels               | mm    | pixels               | mm    |
|                       | 6922                 | 346,1 | 8368                 | 418,4 |

The scope of supply

| Name of component parts of the system                                                                        | Designation of the design document | Quantity, pcs. | Manufacturer/ Note               |
|--------------------------------------------------------------------------------------------------------------|------------------------------------|----------------|----------------------------------|
| Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY:                         |                                    |                |                                  |
| Composition:                                                                                                 |                                    |                |                                  |
| Device for printing digital diagnostic medical images of model MS-460DY                                      |                                    | 1              | Huqiu Imaging (Suzhou) Co., Ltd. |
| Network cable                                                                                                |                                    | 1              |                                  |
| Cable for connecting to a computer network*                                                                  |                                    | 1              |                                  |
| Application instruction**                                                                                    |                                    |                |                                  |
| * Option, available if required.                                                                             |                                    |                |                                  |
| ** The application instruction can be supplied on electronic media. A paper version is available on request. |                                    |                |                                  |

| Name of component parts of the system                                                                        | Designation of the design document | Quantity, pcs. | Manufacturer/ Note               |
|--------------------------------------------------------------------------------------------------------------|------------------------------------|----------------|----------------------------------|
| Device for printing digital diagnostic medical images of models: MS-460DY, MS-760DY:                         |                                    |                |                                  |
| Composition:                                                                                                 |                                    |                |                                  |
| Device for printing digital diagnostic medical images of model MS-760DY                                      |                                    | 1              | Huqiu Imaging (Suzhou) Co., Ltd. |
| Network cable                                                                                                |                                    | 1              |                                  |
| Cable for connecting to a computer network*                                                                  |                                    | 1              |                                  |
| Application instruction**                                                                                    |                                    |                |                                  |
| * Option, available if required.                                                                             |                                    |                |                                  |
| ** The application instruction can be supplied on electronic media. A paper version is available on request. |                                    |                |                                  |

## PRODUCT CLAIMS

---

Any healthcare professional (such as a customer or user) who has a complaint about the equipment, or who is dissatisfied with the performance, longevity, reliability, safety of use, efficiency or performance of this equipment should notify the company Huqiu Imaging (Suzhou)Co., Ltd or its authorized representative in the Russian Federation - Medcorp-MT.

If a serious incident occurs during the use of the MEDICAL DEVICE or as a result of its use, you should report it to the manufacturer and / or its official representative in the Russian Federation at:

Stromynka str., building 19, housing. 2, Э 1 ПИВБ К32 PM3 107076 Moscow, Russia,  
Tel. +7 (495) 969-61-47, [medcorp@medcorp.ru](mailto:medcorp@medcorp.ru)

## INTRODUCTION

---

### Sections:

- INTENDED CONSUMERS AND USERS
- FUNCTIONS AND FEATURES
- EQUIPMENT CLASSIFICATION
- OPERATING RISKS AND PRECAUTIONS
- INSTALLATION
- TRAINING
- COMPLIANCE WITH REGULATIONS AND STANDARDS
- INTERACTION WITH EXTERNAL SYSTEMS
- MARKING
- INSTRUCTIONS FOR USING POWER CABLES
- PATIENT DATA SECURITY
- OPERATING MODES
- LOCAL USER INTERFACE
- TURN ON THE DEVICE FOR PRINTING
- COOLING DEVICE FOR PRINTING
- TURN OFF THE PRINT DEVICE

## **INTENDED CONSUMERS AND USERS**

---

Potential consumers of the Device for printing digital diagnostic medical images of models are patients of medical institutions without age restrictions.

Potential users are qualified personnel of radiographic departments who have passed the appropriate training course.

The term "users" refers to persons who directly work with the equipment, as well as exercise control over its use.

Before operating this equipment, the user must read, understand, take the hint and comply with all warning and instructional labels provided on items of the equipment.

## FUNCTIONS AND FEATURES

---

Device for printing is characterized by the following features and capabilities:

- Thermal printing technology used to obtain high-quality diagnostic images in daylight, has important advantages: no chemicals, washing, drying, lengthy preparations; eliminates the need for a darkened room and the costly disposal of chemical waste. Consumables do not need to be hidden from daylight.

- The compact design of the Device for printing saves space and provides ease of use. The need for maintenance and care is reduced to a minimum.

- Direct thermal printing system provides high quality black and white images: Resolution 320 ppi for MS-460DY, 508 ppi for MS-760DY, 14 bit grayscale for each pixel.

- Films of various sizes can be used.

- The Device for printing input trays are equipped with an RFID reader that automatically scans films and protects the Device for printing if it encounters any media that it cannot recognize.

- Number of input trays. The Device for printing is equipped with two input trays.

- Number of output trays. Device for printing is equipped with one output tray that can be used to handle a variety of film sizes without additional setup.

## EQUIPMENT CLASSIFICATION

|                                                                                                                        |                                                                                                                              |
|------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
| Potential risk class of the MD according to the nomenclature classification of the MD                                  | 1                                                                                                                            |
| IEC 60601-1 type of protection against electric shock                                                                  | Class I equipment - with wired connection                                                                                    |
| IEC 60601-1 degree of protection against electric shock                                                                | The equipment components comply with type B                                                                                  |
| Classification by degree of protection against water penetration and particulates according to GOST R IEC 60601-1-2022 | IPX0 (base is not marked)                                                                                                    |
| Classification by operation mode according to GOST R IEC 60601-1-2022                                                  | Continuous operation                                                                                                         |
| Cleaning                                                                                                               | See Section on cleaning and disinfecting.                                                                                    |
| Disinfection                                                                                                           | See Section on cleaning and disinfecting.                                                                                    |
| Flammable anaesthetic substances                                                                                       | This equipment must not be used in the presence of flammable anaesthetic mixtures with air, oxygen or nitrous oxide          |
| Firmware class according to GOST IEC 62304-2022                                                                        | A                                                                                                                            |
| Waste class according to SanPiN 2.1.3684-21                                                                            | A class                                                                                                                      |
| Russian Classification of Products by Economic Activities 2                                                            | 32.50.50.190                                                                                                                 |
| Climate modification according to GOST R 50444-2020                                                                    | NF category 4.2                                                                                                              |
| Group depending on perceived mechanical influences according to GOST R 50444-2020                                      | Group 2                                                                                                                      |
| MD type code                                                                                                           | 348970                                                                                                                       |
| Software                                                                                                               | HQ_Imager_App_huA24, version 4.44                                                                                            |
| Expected service life                                                                                                  | Five (5) years (subject to regular service and maintenance according to the instructions of Huqiu Imaging (Suzhou) Co., Ltd) |

### Networking capabilities



- The modular design of Device for printing ensures optimal performance in a given network environment.
- Functionally Device for printing fully managed via the network.
- You can control Device for printing using the touch panel located on it or via the network using remote access from a workstation.

### Software Information



- This Device for printing uses the firmware HQ\_Imager\_App\_huA24 Ver. 4.44 developed by Huqiu Imaging (Suzhou)Co., Ltd, China, and is used solely to control the functions of Device for printing.

## OPERATING RISKS AND PRECAUTIONS

|                                                                                                                                                                                                                                                           |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                                                                                                          | CAUTION: Safety when using the printing device is only guaranteed if it has been installed by qualified personnel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|                                                                                                                                                                          | CAUTION: Any medical equipment must be operated by specially trained and qualified personnel.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|                                                                                                                                                                          | CAUTION: Device for printing must be used strictly for its intended purpose and in compliance with the necessary technical requirements. Use of the product for purposes other than those intended or in violation of specifications may result in an emergency situation that could endanger human health, even death (for example, electric shock). In such cases, the manufacturer is not responsible for the consequences.                                                                                                                                                                                                                                                                                                |
|                                                                                                                                                                          | CAUTION: Improper modification, maintenance, or repair of the system can result in personal injury, electric shock, and equipment damage. Safety is only guaranteed if all maintenance and repair activities are carried out by certified field service technicians. Any modification or maintenance of a MEDICAL DEVICE by an uncertified technician is at your own risk and will void your warranty.                                                                                                                                                                                                                                                                                                                        |
|                                                                                                                                                                          | CAUTION: System failure due to hardware or software failure. If this product is used in critical clinical applications, a redundant system must be provided.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|                                                                                                                                                                         | CAUTION: Any images obtained using any printing technology may contain artifacts that can be confused with relevant diagnostic information. If there are doubts about the correctness of the obtained diagnostic information, then for an accurate diagnosis it is necessary to conduct additional studies.                                                                                                                                                                                                                                                                                                                                                                                                                   |
|                                                                                                                                                                        | Clearing a film jam and cleaning the print head can be performed without turning off the power. However, care must be taken and the following instructions must be followed:                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|                                                                                                                                                                        | CAUTION: Observe the following precautions when using and caring for the printing device: <ul style="list-style-type: none"> <li>• Faults in the electrical and mechanical parts must only be eliminated by experienced specialists!</li> <li>• Modification of the equipment is prohibited.</li> <li>• Do not modify or disable the built-in security systems.</li> <li>• Do not block the ventilation openings.</li> <li>• Install the Device for printing so that it can be turned off quickly by unplugging it from the mains.</li> <li>• Before carrying out any maintenance work on the device for printing, always turn off the power to the Device for printing and unplug the power cord from the outlet.</li> </ul> |
|                                                                                                                                                                        | ATTENTION: All warnings, instructions, and safety rules that appear in this document or on parts of the equipment must be strictly observed.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| The product classification according to the IEC 60601-1 standard for medical electrical equipment is intended to be installed outside the patient's immediate vicinity. The parameters of the zone of immediate proximity to the patient are given below. |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |

1. Radius = 1,5 m (EN 60601-1) or 1,83 m (UL 60601-1).
2. Height = 2,5 m (EN 60601-1) or 2,29 m (UL 60601-1).


### Safety Warnings

|  |                                                                                                                  |
|--|------------------------------------------------------------------------------------------------------------------|
|  | CAUTION: Printed images are confidential documents and should only be accessed by authorized personnel.          |
|  | CAUTION: It is a good practice to keep images from being removed from the modal until they are properly printed. |
|  | CAUTION: Images that show film defects are recommended to be reprinted.                                          |

### SYMBOLS AND CONVENTIONS

Symbols are necessary to draw the user's attention to important information, opposite which they are located. Remember the meaning of the symbols. Refer to this section whenever you are unsure of the meaning of symbols or conventions.

| Symbol image | Meaning                                                                                                                        |
|--------------|--------------------------------------------------------------------------------------------------------------------------------|
|              | ATTENTION!                                                                                                                     |
|              | RISK OF ELECTRIC SHOCK!                                                                                                        |
|              | WARNING HOT SURFACE                                                                                                            |
|              | <b>ONLY QUALIFIED PERSONNEL!</b><br>All operations accompanied by this symbol must only be carried out by qualified personnel. |
|              | POWER SUPPLY VOLTAGE TYPE (AC)                                                                                                 |
|              | POWER SUPPLY VOLTAGE TYPE (DC)                                                                                                 |
|              | Grounding                                                                                                                      |
|              | <b>FOLLOW THE INSTRUCTIONS FOR USE.</b><br>FOR CORRECT AND SAFE OPERATION OF THE MEDICAL DEVICE, READ THE INSTRUCTIONS FOR USE |
|              | MANUFACTURER                                                                                                                   |

|                                                                                     |                                                                            |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------|
|    | DATE OF MANUFACTURE                                                        |
|    | SERIAL NUMBER                                                              |
|    | TOP                                                                        |
|    | CE MARKING                                                                 |
|    | IT IS FORBIDDEN TO PUT OTHER LOADS ON THE LOAD IN THE TRANSPORT PACKAGING. |
|    | THE LOAD SHOULD NOT BE ROLLED OR TURNED                                    |
|    | KEEP DRY                                                                   |
|   | TEMPERATURE LIMIT                                                          |
|  | HAZARDOUS WASTE, DO NOT DISPOSING IN HOUSEHOLD WASTE! DISPOSAL REQUIRED    |



## INSTALLATION


Installation and configuration of the Device for printing must be carried out by a local service organization. Some configuration selection tasks (to a limited extent) can be performed by the customer after completing the training. For more information, contact your local service organization.

### Changing the location of the product after installation

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | CAUTION: Be sure to turn it off before moving it.                                                                                                                                                                                                                                                                                                                                                                                     |
|    | CAUTION: All device for printing covers must be closed during transport.                                                                                                                                                                                                                                                                                                                                                              |
|    | CAUTION: Do not lift by the output tray.<br>CAUTION: When moving the Device for printing, consider the stability and structure of the working surface of the table on which the Device for printing will be installed. To avoid interference with ventilation, which may cause overheating. Device for printing should not be placed on soft surfaces. Make sure the Device for printing is placed on a stable table with a hard top. |
|   | CAUTION: Do not move the device from one place to another all the time.                                                                                                                                                                                                                                                                                                                                                               |
|  | CAUTION: Device for printing must be moved by 3 people, or at least 2 people.                                                                                                                                                                                                                                                                                                                                                         |

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                      |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <p>Transporting the Device for printing.</p> <ol style="list-style-type: none"> <li>1. Turn off Device for printing.</li> <li>2. Disconnect the cables.</li> <li>3. Move the Device for printing to a new location (this should be done by 2, preferably 3 people!).</li> <li>4. Connect cables.</li> <li>5. Turn on Device for printing.</li> </ol> |
|  |                                                                                                                                                                                                                                                                                                                                                      |

## SPACE REQUIREMENTS FOR THE ROOM

1. Space required to service the device for printing
2. Space required for normal use

## TRAINING

---

Before starting to work with the MEDICAL DEVICE, the user must receive appropriate training and acquire basic skills for the safe and effective use of the MEDICAL DEVICE. In some countries, the requirements for the training of personnel may have individual specifics.

Users must ensure that they have been trained in accordance with local laws or regulations that are legally binding. Detailed information on personnel training can be obtained from the authorized representative of the manufacturer in the Russian Federation.

The user should pay particular attention to the following information in the system documentation:

- Purpose.
- Intended users.
- Safety instructions.

## LIST OF APPLICABLE STANDARDS

| №  | International standard  | National standard or translation of international standard      | Name of the standard                                                                                                                                                                                          |
|----|-------------------------|-----------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1  | DIN EN 1041:2013        | DIN EN 1041 - 2013                                              | Medical products. Information provided by the manufacturer                                                                                                                                                    |
| 2  | DIN EN ISO 13485: 2016  | GOST ISO 13485-2017                                             | Medical devices - QUALITY MANAGEMENT SYSTEMS - Requirements for regulatory purposes                                                                                                                           |
| 3  | DIN EN ISO 14971: 2019  | GOST ISO 14971-2021                                             | MEDICAL PRODUCTS - Application of risk management to medical devices                                                                                                                                          |
| 4  | DIN EN ISO 15223-1:2016 | GOST R ISO 15223-1-2020                                         | MEDICAL PRODUCTS<br>Symbols used in the labeling of medical devices, on labels and in accompanying documentation<br>Part 1 - Basic Requirements                                                               |
| 5  | IEC 60601-1:2020        | GOST R International Electrotechnical Commission 60601-1-2022   | Medical electrical equipment. Part 1. General requirements for basic safety and essential performance                                                                                                         |
| 6  | DIN EN 60601-1-2:2016   | GOST R International Electrotechnical Commission 60601-1-2-2014 | MEDICAL ELECTRICAL EQUIPMENT - Part 1-2 - General requirements for safety, taking into account essential functional characteristics. Parallel standard. Electromagnetic compatibility. Requirements and tests |
| 7  | DIN EN 60601-1-6:2016   | GOST R International Electrotechnical Commission 60601-1-6-2014 | MEDICAL ELECTRICAL EQUIPMENT - Part 1-6 - General requirements for safety, taking into account essential functional characteristics. Additional standard - USABILITY                                          |
| 8  | DIN EN 62304:2016       | GOST IEC 62304-2022                                             | MEDICAL PRODUCTS - Software. Life cycle processes                                                                                                                                                             |
| 9  | DIN EN 62366-1:2017     | GOST R International Electrotechnical Commission 62366-1-2021   | MEDICAL DEVICES - Part 1 - Design of medical devices for serviceability                                                                                                                                       |
| 10 | —                       | GOST R 50444-2020                                               | Devices and equipment medical. General technical requirements                                                                                                                                                 |

## INTERACTION WITH EXTERNAL SYSTEMS

The Device for printing should only be used in conjunction with equipment or equipment components expressly declared by Huqiu Imaging (Suzhou) Co., Ltd to be compatible with this Device for printing. Modification of the Device for printing is prohibited. To print and obtain clear and dense images in grayscale, the original Thermographic film for medical imaging systems "Dry Film MS-KX410 A/F" with a blue base of various formats 8x10" (20x25cm), 10x12" (25x30cm), 11x14" (28x35cm) and 14x17" (35x43cm), certificate for the medical device No.P3H 2023/21140, manufactured by Huqiu Imaging (Suzhou) Co., Ltd should be used in conjunction with the Device for printing digital diagnostic medical images of models: MS- 460DY, MS-760DY, manufactured by Huqiu Imaging (Suzhou) Co., Ltd. Any other types of film from any other manufacturers are prohibited due to lack of compatibility.

This Device for printing is a desktop digital thermal device for printing digital diagnostic medical images. This means that this Device for printing can be connected directly to an (existing) local Ethernet network without the need for additional, incl. auxiliary equipment. Also, Device for printing is fully compatible with the DICOM protocol. In this regard, the standard DICOM protocol can be used as a network protocol: thus, Device for printing is able to process and print tasks within the DICOM protocol without using additional, incl. auxiliary equipment.

### Safety Compliance

Peripheral equipment connected through any interfaces must be certified according to the applicable IEC regulations (eg IEC 62368 for data processing equipment or IEC 60601-1 for medical equipment). In addition, all systems must meet the requirements for medical equipment (ME) in accordance with IEC 60601-1. Individuals or firms that connect additional equipment to an input signal node or an output signal node are setting up a MEDICAL DEVICE and are therefore responsible for ensuring that the system meets the IEC 60601-1 medical equipment (ME) requirements. If you have any questions, please contact the authorized representative of the manufacturer in the Russian Federation.


## MARKING

### Marking of the MEDICAL DEVICE

The Device for printing digital diagnostic medical images is labeled with the following information:

- Logo of the manufacturer;
- Device name;
- Model;
- Characteristics of the device (electrical)
- Serial number;
- Issue date (year, month);
- Manufacturer, place of production.

Marking location:

### Packaging marking

The symbols on the boxes are:

|  |                                                                            |
|--|----------------------------------------------------------------------------|
|  | TOP                                                                        |
|  | IT IS FORBIDDEN TO PUT OTHER LOADS ON THE LOAD IN THE TRANSPORT PACKAGING. |
|  | THE LOAD SHOULD NOT BE ROLLED OR TURNED                                    |
|  | KEEP DRY                                                                   |

## **INSTRUCTIONS FOR USING POWER CABLES**

---

- When operating the device, use only the supplied power cable set.
- Do not use the supplied power cable set with other electrical appliances.

## PATIENT DATA SECURITY

The user must ensure that the legal rights of patients are respected, as well as take measures in relation to ensuring the appropriate security of patient data.

The user defines who has access to patient data in certain situations.

The user must define a strategy for dealing with patient data in emergency situations.

### Data storage

The user must ensure that confidential patient information cannot be retrieved from the built-in storage media (magnetic hard drive, solid state drive, etc.) before the Device for printing (or built-in storage device) leaves the office.

### Operating environment requirements

These information security and privacy operating environment requirements established in accordance with paragraphs 17(4) and 18(8) of Annex I to the EU Medical Devices Regulation 2017/745 must be implemented and applied by the Client (User) in connection with the use of a MEDICAL DEVICE by Huqiu Imaging (Suzhou) Co., Ltd. These minimum requirements are intended to prevent unauthorized access that could prevent the device from functioning as intended.

- Perimetric firewalls should be provided and properly configured to ensure that communication between medical devices and external resources is guaranteed to be prohibited or restricted to only those data absolutely necessary for the proper functioning of medical devices.
- Network Intrusion Detection/Prevention Systems (NIDS/NIPS) perimetric detection/prevention systems should be provided and properly configured to provide successful early warning of attempted attacks or threats to medical devices and to prevent threats to medical devices.
- Medical device configurations must specify a Network Time Protocol (NTP) server in order to synchronize the time in the audit logs with the NTP server time.
- Medical devices must belong to an isolated network segment in which communication between medical devices and systems is limited to only the data required for the devices to function.
- Internal firewalls should be used to improve network segmentation and further restrict communications between medical devices and the systems (internal and external) with which they must communicate.
- You should back up your medical device configurations to a separate secure device.
- Security controls should be arranged to ensure that physical access to medical devices is restricted to authorized individuals, and to ensure that devices are not physically stolen.
- An incident response plan should be developed that details the responsibility and response in the event of an incident. Personnel involved in the implementation of the incident response plan should be trained to ensure proper and effective response.
- A formal user registration and deregistration process should be arranged to ensure proper management of access rights to medical devices.

- Users should be assigned unique accounts for medical devices.
- User access rights to medical devices should be reviewed as necessary for appropriateness and adjusted at regular intervals, but at least annually.

## OPERATING MODES

---

Device for printing can operate in two modes: operator mode, service mode

### **Control modes (local and remote)**

The printing device can be controlled locally from the touch panel of the printing device, and via the network using remote access from the workstation.

### **Operator mode**

The operator mode includes all basic functions that are intended for radiologists without special technical skills:

- Creation of pictures suitable for diagnostics;
- Film loading;
- Ensuring the normal functioning of the printing device.

All operator mode functions are described in the operating instructions.

Access is possible using both a local touch panel and a connected remote PC (with password protection).

### **Service mode**

The service mode functions are intended for use by service technicians. Access to the service mode is password protected.


Access is possible using a connected remote PC.

## LOCAL USER INTERFACE

Device for printing interacts with the user through the following controls:

- button "Power/Reset";
- control display;
- light signals.
- sound signals.

1. Film output tray
2. Touch display control
3. Power/Reset Button
4. Pause button
5. Film input trays

### CAUTION:

Never open Device for printing while printing. Always follow the instructions on the display!

### Messages

Under certain circumstances, the display will show a message, see figure below, warning of an error. Such messages inform the user about the occurrence of errors or the impossibility of performing the requested action/operation.

**hu.q**

Error Log

The user should read these messages carefully. They contain information on measures/next steps to be taken in a given situation. Such measures / actions are the implementation of certain operations that eliminate the problem that has arisen, or contacting the local service organization. Detailed information on the contents of the messages can be found in the service documentation provided to service personnel.


*Note:* A list of error codes is provided in Appendix A of the operating instructions.

#### Control buttons

|              |                                                                                                                                                                       |              |                                                                                                                                   |
|--------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------|-----------------------------------------------------------------------------------------------------------------------------------|
| <b>POWER</b> | Power                                                                                                                                                                 | Pause button | <ul style="list-style-type: none"><li>• Turns on/off the printing device.</li><li>• Stops the print queue of a picture.</li></ul> |
| <b>PAUSE</b> |                                                                                                                                                                       |              |                                                                                                                                   |
|              | <b>CAUTION:</b><br>DO NOT press the Power button before completing the procedure for stopping the printing process while the Device for printing is printing on film. |              |                                                                                                                                   |

#### Back panel

|                                                                                    |  |                                                                                                                                                                                                                                                                                                                                                                      |
|------------------------------------------------------------------------------------|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  |  | The back of the Device for printing has one slot and three connectors:<br>1. USB port socket<br>For inserting a USB flash drive when you need to install software, make a backup, etc.<br>2. Network connector<br>Serves for connection to a computer network of medical institution.<br>3. Power connector<br>Used to connect the power cable<br>1.4. Power button. |
|------------------------------------------------------------------------------------|--|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|

#### USB Inserting a USB flash drive

Device for printing has a USB port on the back.

The USB flash drive must be formatted using the FAT32 file system. When connecting a USB flash drive with a different file system, an error message is displayed.

#### Light and sound signals

Device for printing presents information about its status using light and sound signals on the touch control panel and sound. By the color and sound of the signal, you can determine the nature of the response of the device:

- Green - Ready to print, no film in the tray.
- Yellow – Film is in the tray, the thermal head is heating or printing is in progress.
- Red - error (for example, the film is not loaded correctly)

#### Status indicator




| Color / Signal |                      | Condition             | Action                     |
|----------------|----------------------|-----------------------|----------------------------|
| Green          | Lights up constantly | Ready to go (waiting) | Keep going.                |
| Red            | Flashing             | Warning               | Check messages on display. |
|                | Lights up constantly | Error                 |                            |

The Device for printing presents information about its status with beeps. By the duration of the sound signal, you can determine the nature of the system's response to commands from the keypad.

- A **short** beep indicates that you have received a command when you press a touch key and have started executing it, for example, when you turn on the Device for printing.
- A **long** beep indicates that a printing error or other error has occurred, an error code is displayed on the touch control panel.

## Display

|                                                                                                                         |                                                                                                                                                                                                                                                                                           |
|-------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                        | <p>Touch control display, used to control image printing and to enter numeric or alphanumeric data. The figure below shows the display layout used in the INSTRUCTIONS FOR USE.</p>                                                                                                       |
|                                      |                                                                                                                                                                                                                                                                                           |
| <p><b>Main menu</b></p>              | <p>The touch screen provides text output. The main functions of the display are shown below.</p>                                                                                                                                                                                          |
| <p><b>Print task management</b></p>  | <p>During printing, this screen displays current print tasks. To go to the confirmation window for deleting a task, click on the respective task and then click «OK» to delete the task or «Cancel» to return to the print task window. Note: The first print task cannot be deleted.</p> |

|                                                                                                                                                                                                                                                                                                                                                      |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                      |     |     |     |     |       |     |     |     |     |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----|-----|-----|-----|-------|-----|-----|-----|-----|
| <p><b>Sample image printing</b></p> <p><b>hu.q</b></p> <p>Print Size</p> <div></div> <div><div>8 * 10</div><div>10 * 12</div><div>11 * 14</div><div>14 * 17</div></div>            | <p>Selecting a picture for Sample image: Select a picture for Sample image. (Note: If there is no film in the corresponding tray, the following window will not appear on the display.)</p> <p><b>hu.q</b></p> <p>Picture</p> <div></div> <table><tr><td>CT1</td><td>CT2</td><td>CT3</td></tr><tr><td>CT4</td><td>X-Ray</td><td>MR1</td></tr><tr><td>MR2</td><td>MR3</td><td>MR4</td></tr></table> | CT1 | CT2 | CT3 | CT4 | X-Ray | MR1 | MR2 | MR3 | MR4 |
| CT1                                                                                                                                                                                                                                                                                                                                                  | CT2                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | CT3 |     |     |     |       |     |     |     |     |
| CT4                                                                                                                                                                                                                                                                                                                                                  | X-Ray                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                | MR1 |     |     |     |       |     |     |     |     |
| MR2                                                                                                                                                                                                                                                                                                                                                  | MR3                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  | MR4 |     |     |     |       |     |     |     |     |
| <p><b>Service mode</b></p> <p><b>hu.q</b></p> <p>Service</p> <div></div> <div><div>Sensor Status1</div><div>Sensor Status2</div><div>Motor Control</div><div>Error Log</div></div> | <p>In the service mode, the operator can control the functioning of the sensors responsible for the positioning and placement of the film in the trays, the degree of pressing of the thermal head of the printing device and its temperature, control the drive of the print head, the drives of the film clamp and transport. When an error occurs, a special message with an error code is displayed on the display.</p>                                                                                                                                          |     |     |     |     |       |     |     |     |     |
| <p><b>System Information</b></p> <p><b>hu.q</b></p> <p>Sys Info</p> <div></div> <div><div>SN:</div><div>PC Ver:</div><div>Soft Ver:</div><div>AEtitle:</div><div>IP:</div><div>Port:</div></div>                                                                  | <p>This display shows system information such as software version, IP address</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |     |     |     |     |       |     |     |     |     |

#### Data input

When entering numeric or alphanumeric data, always adhere to the following guidelines:

- Only alphanumeric data can be entered.
- During input, the data is displayed in inverse mode.
- You can increase the alphanumeric value you enter using the «**TouchPad Up**» key. Going from 9 to 0 of the value of one number will also increase by one the adjacent number on the left, taking into account the allowable range.
- You can decrease the alphanumeric value you enter using the «**TouchPad Down**» key. Going from 0 to 9 of the value of one number will also decrease by one the adjacent number on the left, taking into account the allowable range.
- To repeat a move command multiple times, you can press and hold the corresponding key.
- To enter a command in a menu, press the «**Confirm**» key.
- Command entry is confirmed by a short beep.
- Device for printing emits a long beep if a key has been pressed that is currently inactive.

## TURN ON THE DEVICE FOR PRINTING

*Note:* Read the safety instructions before turning on the Device for printing.

The following describes how to properly turn on and check the operation of the Device for printing.

1. Check if the power cable is connected and turn on by pressing the **Power** key on the back of the Device for printing.


2. After connecting the device for printing to the network, start by pressing the «**POWER**» button on the front of the device for printing control touch display.

The start screen appears on the touch screen.


After starting, the following message is displayed. For some time, the progress bar will show the passing of the self-test. Next, a screen will appear with the following information:

It takes approximately 11 minutes to prepare Device for printing for use in print mode. Approximately 6 minutes later, a message indicating that Device for printing is ready to print appears, after which you can send print tasks, but it will take another five minutes before the head warms up.

If you send print tasks during the specified five-minute period, then during the remaining time it will perform calculations as part of processing the received task, after which a message will appear on the display stating that it is warming up.

3. Device for printing is ready to go:

If the front panel display shows «**READY**», the status indicator is permanent green.

4. Ensure that the correct supplies are loaded.

*Note:* If the task status indicator contains a warning or an error, see the troubleshooting section.

## COOLING DEVICE FOR PRINTING

---

If, due to certain circumstances (for example, prolonged operation in print mode), there is a significant increase in the temperature of Device for printing, it automatically switches to cooling mode. A message will appear on the display while Device for printing is in cooling mode.

## TURN OFF THE PRINT DEVICE

When turning off, remember that before you complete the task, it will finish the current print task. Print tasks that are in the queue will be stored internally and will be resumed the next time Device for printing is turned on.

Go from the main menu of Device for printing by the corresponding link «Close Imager»

hu.q

Menu

Mission Manage

Sample Print

Service Mode

Close Imager

hu.q

Close Imager

OK

Cancel

To turn off press OK.

*Note:* After turning off, the device enters standby mode. To disconnect the device from the power source, first unplug it using the power switch and unplug the power cord from the power outlet.

*Note:* Safe shutdown is guaranteed when a device for printing completes a current print task. After completion an operator must act in accordance with a section TURN OFF THE PRINT DEVICE and then disconnect the device from the power supply.

## BASIC OPERATIONS (USER MODE)

This section describes the basic principles of working with Device for printing. After completing this section, the operator will be able to print diagnostic images. No special technical skills are required.

All major operator mode functions can be directly accessed by pressing a single key on the keypad.

| Function / task                     | Description                                                                                                                                        |
|-------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------|
| PAUSE PRINT TASK                    | Used to pause print tasks that are in the queue. Device for printing will end the current print task, but will not proceed to the next print task. |
| VIEW THE NUMBER OF FILMS            | Used to determine the number films left in the trays.                                                                                              |
| VIEW THE DEVICE FOR PRINTING STATUS | Displays information about current processes/events taking place in Device for printing                                                            |

### Sections:

- PRINT TASK MANAGEMENT
- PAUSE PRINT TASK
- VIEW THE NUMBER OF FILMS
- DELETE PRINTING TASK
- PRINT SIZE
- FILM LOADING

## PRINT TASK MANAGEMENT

- Checking the PRINT TASK
- PRINT TASK Local Screen

### Checking the PRINT TASK


The user can always check the status of print tasks in remote mode from the operator's menu via the network using remote access from a workstation.

**Note:** Keep in mind that multiple films may be required for a single task. Depending on the imaging modality being used and the actual settings, several films can be grouped into one folder and sent to print as one single task. For more information on the modalities used to acquire images, see the INSTRUCTIONS FOR USE.

Print tasks sent to the network are processed in the order in which they are received. New tasks put in the queue receive the status of "waiting".

As soon as the last film of any task is ejected from the output tray, the next task for which calculations have been performed starts printing.


1) Number of current print tasks

### PRINT TASK Local Screen

During printing, the local display shows the "printing" icon:

The screen displays information about the task that is currently being printed:


During printing, this screen displays current print tasks.


## PAUSE PRINT TASK

During the printing process, the user can pause the print task processing at any time by entering the operator mode menu.

While printing is in progress, the «PRINT» icon is displayed:

To go to the confirmation window for deleting a task, click on the respective task and then click «OK» to delete the task or «CANCEL» to return to the print task window.

## VIEW THE NUMBER OF FILMS

During the printing process, the user can at any time view information on the number of films remaining in the trays, which is displayed on the display of the device.


### CAUTION:


While printing, do not open the input trays to view the number of films loaded, but instead follow the instructions below:

To view the number of films in the trays, follow these steps:


While printing is in progress, the «PRINT» icon is displayed:

Pause printing and press the «BACK» key to enter the operator menu. The «OPERATOR MENU» will be opened on the display:

If it is in the process of printing, it will first complete the current print task. Information about the number of film sheets is in the size line of the film used when printing the image:

- 1) Film size
- 2) Number of films

*Note:* If a print task cannot be completed because there is no film in the tray being used, the Device for printing checks if there is any film in the second tray and that it can be used for the current and other pending print tasks. Device for printing will skip print tasks that it cannot complete at the moment and return to doing them later.

## DELETE PRINTING TASK

---


*Note:* The user can only delete printing tasks from the task queue using remote access from a workstation (in remote operator mode).

## PRINT SIZE

The MEDICAL DEVICE is compatible with thermal imaging films for medical imaging systems in four formats (8x10 inches (20x25 cm), 10x12 inches (25x30 cm), 11x14 inches (28x35 cm), 14x17 inches (35x43 cm), three of 8x10" (20x25 cm), 10x12" (25x30 cm), 11x14" (28x35 cm) are placed in the top tray. 14x17" (35x43 cm), 11x14" (28x35 cm) films are placed individually in the bottom tray. Each section has two cartridge stops. When printing on 8x10" (20x25 cm), 10x12" (25x30 cm), 11x14" (28x35 cm) films, the cartridge stops must be secured in positions A1 and A2 as shown in the figure below.

When printing on 14 x 17 inch (35 x 43 cm) film, the cartridge stops must be locked in the B1 and B2 positions. When installing cartridge stops, the round latch must be fixed on the left stop in the left slot of the cartridge drawer and then moved to the left to lock it. The same process should be repeated on the right side.

### The location of the film in the trays.

## FILM LOADING

This section describes how to load film with the appropriate specifications.

Into Device for printing the following Film can be loaded: 8x10 inches (20x25 cm), 10x12 inches (25x30 cm), 11x14 inches (28x35 cm), 14x17 inches (35x43 cm).

*Note:* To change the film size, you must change the tray configuration.

Device for printing reports that the film has ended in several ways:

- The status indicator on the touch panel flashes (red),
- The display shows a message that the input tray is empty.

The film loading procedures for the top and bottom trays are the same. In the examples below, it is assumed that film needs to be loaded into the bottom input tray.


**CAUTION:**

Do not add additional sheets of film to the film cartridge currently in use. The empty input tray can only load a full film cartridge.

### Film Loading Procedure


1. Open the empty input tray.

**CAUTION:**

To avoid film jams, the input tray must be fully opened.

2. Take the film cartridge, sealed in a metallized bag, and open it.


*Note:* Make sure that the type of film in the cartridge matches the label on the tray! If you are using a different type of film, we recommend changing the label on the tray.


*Note:* All actions with the package are easier to do by laying it on the table. Before doing this, make sure that there is no dust on the table!

3. Partially remove the metallized bag.

4. Now the film cartridge can be placed in the input tray and the metallized bag can be completely removed.

5. Remove the white protective sheet.  
6. Close the input tray.

Document version v2.20.09.23


|                                                                                   |                                                                                                                                                                                                                                                                                                                                                                                 |
|-----------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|  | <i>Note:</i> Device for printing continues printing as soon as the tray is closed.                                                                                                                                                                                                                                                                                              |
|  | <i>Note:</i> If there is film in the input tray, you cannot load any other size film in the input tray. Changing the film format during printing increases the risk of dust ingress, which can damage the thermal head.                                                                                                                                                         |
|  | <i>Note:</i> If a print task cannot be completed because there is no film in the tray being used, it will check for the presence of film in the second tray and that it can be used for the current and other print tasks in the queue, and skips print tasks that it cannot complete in the current moment, and will return to their implementation later.                     |
|  | <i>Note:</i> If the film cartridge is not correctly loaded into the input tray, it will print the image on the reverse side of the sheet (which does not have an emulsion coating). This will cause a film jam. The operator will be prompted to print the image on the back of the film sheet and prompted to clear the film jam and verify that the film is loaded correctly. |
|  | <b>CAUTION:</b><br>The user should not reuse film that is caught in a jam.                                                                                                                                                                                                                                                                                                      |

## PROPER POSITIONING OF THE FILM IN THE LOADING TRAY


You can check the correct positioning of the film sheets in the tray by looking at the near-left corner of the film. The rounding of this corner has a recess with a notch.

It is this corner that must **ALWAYS** be in the near left corner of the tray when loading the film cartridge into the tray of the Device for printing (with the imaging surface at the bottom).


When the input tray is closed, it reads the film IDs and automatically makes the necessary settings. The identification token is located on the back of the film cartridge. The figure below shows the film cartridge on the back.

## ENU STRUCTURE

|                    |                                                                                                                                                                                |
|--------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Main menu          | Main Menu                                                                                                                                                                      |
| Mission Manage     | Viewing and deleting a print task                                                                                                                                              |
| Sample Print       | Printing Test Images                                                                                                                                                           |
| Service Mode       | Service menu:<br>View the status of sensors of Device for printing, error log, manual control of engines Device for printing. (Used when performing technical service)         |
| Close Imager       | Turn off the Device for printing                                                                                                                                               |
| 8*10<br>10*12      | Sample size selection: Select the sample size to be printed                                                                                                                    |
| CT1<br>CT2<br>MR2  | Sample image selection: Select a sample image for printing (Note: If there is no film in the corresponding field, you will not be able to enter the image selection interface) |
| Sensor Status1     | View sensor status1                                                                                                                                                            |
| Sensor Status2     | View sensor status2                                                                                                                                                            |
| Motor Control      | Manual Calibration of Engines Device for printing                                                                                                                              |
| Error Log          | Error Status of Device for printing                                                                                                                                            |
| System information | System Information: serial number of Device for printing, software version, network settings.                                                                                  |

Document version v2.20.09.23

## QUALITY CONTROL

Regular quality control is essential to determine and maintain consistent image quality.

The Device for printing has an automatic quality control procedure that meets the requirements for gray scale reproducibility tests, as defined by the international standard IEC 1223-2-4.

Local codes and regulations may require the use of other methods.

The quality control procedure consists of quality control of the test image.

- Before starting work, it is necessary to print a test image, which will be used later to control and check the image quality (annually).



*Note:* Follow the print quality check steps for each input tray that uses thermographic film for medical imaging systems.

The results of the checks are recorded in the quality control tables (for example, the Geometry Control Table, see the appropriate section.)

## BEST IMAGE FOR QUALITY CONTROL

SMPTE Image Quality Control Test Example, the printed image must meet the following specifications:

A-width  
B-length

(1) Fields 1 (4X) contain 6 squares with horizontal and vertical black and white stripes (width 1.2 and 3 pixels)

Check: Each of the 6 squares should show a stable distance between the black and white stripes (interference must be eliminated, as this leads to an uneven distance between the stripes)

(2) Fields 2 (4X) contain the same pattern, but with less contrast.

Check: All 6 squares must be visible.

(3) Fields 3 contain a certain number of squares with the same density.

Check: All individual squares must have a constant uniform density, without stripes and dots.

(4) Field 4 contains a black square at 100% optical density with a smaller square inside having 95% density.

The field contains an area with two different densities, i.e. a small square inside a large one, with a slight difference in blackening density, should be visible.

(5) Field 5 contains a light square (0% density) with a small gray square inside (5%). The field contains an area with two different densities, i.e. a small square inside a large one, with a slight difference in blackening density, should be visible.

(6) Field 6 contains a black square with a maximum density of 100%

Check: Square should show Dmax (default is 3.2 O.D.)

(7) Fields 7 contain light squares with a density of 0%

Check: The square should indicate the degree of blackening of the veil (0.2 O.D. for blue-based film)

### CAUTION:

Good viewing conditions are essential for the correct interpretation of both diagnostic and test images. Set the light intensity (luminosity) to between 2000 and 4000 cd/m<sup>2</sup> (color temperature between 4500 and 6500°K). Use a magnifying glass and shutters to collimate the image. The ambient light level should be low.

## CHECKING THE ALLOWABLE SPATIAL RESOLUTION, THE LEVEL OF INTERFERENCE AND THE VISIBILITY OF DETAILS WITH LOW CONTRAST

To check the acceptable spatial resolution, interference level, and visibility of low contrast details, do the following:

1. Print a quality control test image, or use a previously printed image that you use in your daily work.
2. Visually inspect the QC test image for interference: there should be no interference that significantly interferes with visual perception.
3. Check the visibility of low contrast details at both the upper (100 / 95%) and lower (0 / 5%) limits of the density range. You should see field 5, which contains a light square (0% density) with a small gray square inside (5%). The field contains an area with two different densities, i.e. a small square inside a large one, with a slight difference in blackening density, should be visible. Field 4 should have a black square with 100% optical density with a smaller square inside having 95% density. The field contains an area with two different densities, i.e. a small square inside a large one, with a slight difference in blackening density, should be visible.

### CAUTION:

In the event of interference, insufficient spatial resolution, or errors in any other recommended quality control check, the cause of the problem must be found and corrected before any diagnostic images can be printed on the Device for printing.

## DETERMINATION OF THE BASIC VALUES OF THE GEOMETRIC PARAMETERS OF THE IMAGE

To determine the basic values of the geometric parameters of the image, follow the instructions below:

1. Print a test image for quality control or use a previously printed one.
2. To determine the basic values of geometric dimensions, measure the sides A and B of the square on the test image, see figure.

### CAUTION:

Side A must be measured from the left edge of the left side to the right edge of the right side, and side B - from the top edge of the top side to the bottom edge of the bottom side.

We strongly recommend using a 30 cm ruler with 0.5 mm graduations for this.

3. Record these values as base dimensions  $A_{ref}$  and  $B_{ref}$  in the table below ("Table of geometric parameters control").

These tables are used for annual image quality checks.

4. Save this film for future use.

### Table of geometric parameters control

Test frequency:

Annually or as needed

Serial number \_\_\_\_\_

Input tray:

| Initial dimensions<br>Date: |  | Measured dimensions<br>Date: |  | Comparability |  | Aspect Ratio |  |
|-----------------------------|--|------------------------------|--|---------------|--|--------------|--|
| $A_{ref}$                   |  | A                            |  | $A / A_{ref}$ |  | A/B          |  |
| $B_{ref}$                   |  | B                            |  | $B / B_{ref}$ |  |              |  |

| Initial dimensions<br>Date: |  | Measured dimensions<br>Date: |  | Comparability |  | Aspect Ratio |  |
|-----------------------------|--|------------------------------|--|---------------|--|--------------|--|
| $A_{ref}$                   |  | A                            |  | $A / A_{ref}$ |  | A/B          |  |
| $B_{ref}$                   |  | B                            |  | $B / B_{ref}$ |  |              |  |

To detect deviations in image sizes and aspect ratios, it is necessary to carry out this check once a year.

1. Do a daily check first.
2. Take a QC test image as part of the weekly inspection and measure sides A and B of the square.
3. Record the dimensions of sides A and B in the table.
4. Compare the measured values A and B with the base values of the  $A_{ref}$  and  $B_{ref}$  parameters from the Table

The difference between the measured values of sides A and B and the basic values  $A_{ref}$  and  $B_{ref}$  should not exceed 1.0%.

5. Examine the image for distortion.

6. Calculate the aspect ratio by dividing A by B.

The result should be equal to  $1 \pm 0.01$

### CAUTION:

If the image size or shape distortion is out of range, contact your service organization to resolve the problem.

## SCHEDULE OF PREVENTIVE CARE ACTIVITIES

The design of the Device for printing unit ensures reliable and trouble-free operation. Care and cleaning require minimal effort from the user.

The following pages describe procedures for cleaning the MEDICAL DEVICE.

| Periodicity                                             | What should be done?        |
|---------------------------------------------------------|-----------------------------|
| Episodically                                            | "Cleaning and Disinfection" |
| When there are signs of deterioration in image quality. | "Cleaning the Print Head"   |

Cleaning the print head should only be performed when there is deterioration in print quality. For questions regarding service procedures, always contact your local service representative.

## SAFETY INSTRUCTIONS

### CAUTION:

To prevent damage to the Device for printing, observe the following precautions:

- Do not lubricate.
- Do not disassemble.
- Do not touch the print head fuser array.
- Always turn off and unplug the power cord before doing any work inside the Device for printing.

*Note:* Clearing a film jam and cleaning the print head can be performed without turning off the power.

## CLEANING AND DISINFECTION

To avoid contamination of personnel, patients and equipment, all instructions and procedures must be strictly followed. Take general precautions against possible contamination of the Device for printing. See below for more information on cleaning the device.


To clean the outside of the Device for printing:

1. Turn off.
2. Unplug the power cord from the electrical outlet.
3. Wipe the outside of the case (body) with a clean, soft, damp cloth.

Use a lint-free cloth and special products for monitor screens. Be careful not to scratch the touch control panel.

Cleaning and disinfection of the installation should be carried out according to MU-287-113 with any disinfectant solution that does not contain water and is approved for use on the territory of the Russian Federation (for example: 0.5% alcohol solution of chlorhexidine).

Do not use cleaning materials with a high alcohol content. Do not use corrosive detergents, abrasive detergents or solvents. Use only explosion-proof disinfectants. Never use ammonia-based cleaners.

|                                                                                     |                                                                                                                                                                                                                                               |
|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <b>CAUTION:</b><br>Do not allow liquids to enter the Device for printing                                                                                                                                                                      |
|   | <b>ATTENTION!</b><br>The equipment must not be used in the presence of flammable anesthetics and/or the presence of flammable disinfectants. When using flammable disinfectants, wait for the combustible gases to dissipate, do not turn on. |
|  | <i>Note:</i> To clean, do not open the Device for printing case. The internal parts of the print device are not intended to be cleaned by the user.                                                                                           |

4. Connect to the mains and turn it on.

## **REQUIREMENTS FOR THE TECHNICAL SERVICE AND REPAIR OF A MEDICAL DEVICE**

If the equipment fails or if the equipment does not work properly, you must contact the service department.

Repair of the printing device at the consumer's place is carried out by the manufacturer's service department or an organization authorized by the manufacturer and licensed to carry out this type of work.

Current repairs include work on diagnostics, elimination of identified defects and repair of the printing device with the replacement of faulty components, parts, assemblies, verification in operation, which are carried out at the request of the Consumer.

The Manufacturer, upon request of the service organization, will provide wiring diagrams, component specifications, calibration instructions and other information required by Authorized Personnel.

Device for printing needs regular checks and service. The recommendations below are aimed at improving the reliability and safety of equipment operation.

Technical service includes checks and inspections that must be carried out by authorized personnel.

It is necessary to train the personnel working on the printing device to perform daily checks. Technical service must only be carried out by manufacturer authorized technicians.

The calendar average service life of the device must be at least 5 years. The limit state criterion is the impossibility or economic inexpediency of restoring the device by repair.

## **CHECKING AND CONTROL OF THE SYSTEM BY PERSONNEL USING THE DEVICE**

### **Carrying out daily maintenance before using the equipment**

Perform a visual inspection of the technical condition of the equipment (checking the presence of screens, fences, protective devices, warning alarms, etc.).

Check visually the connection to the grounding equipment.

When turning on the equipment, check that the Device for printing has finished booting without any error or failure messages.

Visually check the readiness of the Device for printing for use, the initial positions of the controls, the condition of the indication and signaling devices.

During the operation of the Device for printing, visually check the absence of violations in the operation of the hardware components of the Device for printing, error messages or failures.

Visually check for contamination, liquid organic residues.

Cleaning and disinfection is carried out within the framework of the relevant section of these INSTRUCTIONS FOR USE.

## TECHNICAL SERVICE OF THE SYSTEM AT THE CONSUMER

Technical Service of the Device for printing at the consumer is performed by the manufacturer's service department or an organization authorized by the manufacturer and licensed to carry out this type of work.

Technical Service is carried out to ensure reliable operation and technically sound condition during the entire period of operation of the Device for printing after the end of the warranty period.

Technical Service of the Device for printing consists of:

- checking the technical condition, its parts, components and assemblies;
- cleaning of the thermal head;
- tightening loose fasteners;
- checking the condition of grounding units, cables, connecting conductors, other switching devices;
- verification of the control and management bodies for integrity;
- instrumental control and adjustment of the main technical characteristics:

Cleaning the print head should only be performed when there is deterioration in print quality. Generally, the print head should be cleaned after 100-200 films have been printed. To clean the print head, use a non-woven cloth soaked in water or absolute alcohol.

- image quality control

The Device for printing has an automatic quality control procedure that meets the requirements for gray scale reproducibility tests, as defined by the international standard IEC 1223-2-4. Local codes and regulations may require the use of other methods.

The quality control procedure consists of two main steps:

- Before starting work, it is necessary to carry out a printout of the SMPTE control test and check the quality of the image in accordance with the instruction manual.
- After the test, you need to make sure that there is no interference and that details with low contrast are visible.

## **OUNTING AND INSTALLATION REQUIREMENTS**

---

Mounting, commissioning, calibration and dismantling of the Device for printing must be carried out by authorized specialists only.

When putting the Device for printing into operation, the personnel of the medical institution undergoes a detailed briefing by the manufacturer's specialists or its authorized specialists.

Timely technical service should only be carried out by qualified service personnel.

With proper technical service, it remains operational throughout its entire service life. The performance of the Device for printing is confirmed by periodic checks.

Device for printing is intended for use exclusively by expert users, i.e. trained professionals. Operating personnel should regularly attend refresher courses to work on this type of equipment.

Space requirements for the maintenance and operation room are specified in the relevant section of this INSTRUCTIONS FOR USE. When installing or moving the Device for printing, the appropriate section of this INSTRUCTIONS FOR USE must be observed.

## THE PROCEDURE AND CONDITIONS FOR THE RECYCLING OR DISPOSAL OF THE MEDICAL DEVICE

---

At the end of the full service life, Device for printing must be decommissioned, dismantled and disposed of by an organization licensed to carry out these works.

Decommissioning provides for a set of administrative and technical actions aimed at ensuring the safety of personnel, the public and environmental protection at all stages of operation. All existing laws and regulations must be complied with and disposal work must be carried out by a certified contractor. Be sure to delete all data stored on discs or other media prior to disposal so as not to disclose private information.

Collection, use, disposal, placement, storage, transportation, accounting and disposal of medical waste must be carried out in compliance with the requirements of the Sanitary Rules SanPiN 2.1.3684-21 section X. Device for printing is disposed of as epidemiological safe waste, similar in composition to solid household waste - **class A**.

## **REQUIREMENTS FOR ENVIRONMENTAL PROTECTION WHEN USING THE MEDICAL DEVICE**

This MEDICAL DEVICE does not harm people and the environment when used for its intended purpose, transported, stored and disposed of in accordance with the requirements of the regulatory, technical and operational documentation of the manufacturer.

The MEDICAL DEVICE is accompanied by the  WEEE symbol. The aim of the Waste Electrical and Electronic Equipment (WEEE) directive is to reduce the accumulation of electrical and electronic waste through recycling and other forms of reuse. According to regulations, it is necessary to ensure the collection, processing and reuse of such types of waste. Due to the peculiarities of national laws, different requirements may apply in different countries. The WEEE mark on the product and accompanying documents prohibits the disposal of end-of-life electrical and electronic products with household waste. For details on how to return this product for recycling, please contact your local service provider and/or authorized representative. The use of recycled materials helps conserve natural resources.

## **WARRANTY OBLIGATIONS AND RECLAMATION OF THE MEDICAL DEVICE**

---

Any healthcare professional (such as a customer or user) who has a complaint about the equipment, or who is dissatisfied with the performance, longevity, reliability, safety of use, efficiency, or performance of this equipment should notify the company Huqiu Imaging (Suzhou)Co., Ltd. If the malfunction of the Device for printing has caused or contributed to serious harm to the patient's health, Huqiu Imaging (Suzhou)Co., Ltd should be informed immediately by phone, fax or send a notification by mail to the address below:

- Huqiu Imaging (Suzhou)Co., Ltd
- Medcorp-MT, Stromynka str., building 19, housing. 2, Э 1 ПИВБ K32 PM3 107076 Moscow, Russia  
tel. +7 (495) 969-61-47, [medcorp@medcorp.ru](mailto:medcorp@medcorp.ru)

Support and Service Huqiu Imaging (Suzhou)Co., Ltd. - addresses and phone numbers of local support and service representatives are given on the Internet on the website of the authorized representative in the Russian Federation of Medcorp-MT [www.medcorp.ru](http://www.medcorp.ru)

The manufacturer's warranty is 12 months.

The approximate life of the product with appropriate regular inspection and technical service is 5 years.

Warranty shelf life of MD in a transport package is up to 1 year.

## HIGH-FREQUENCY RADIATION AND PROTECTION INFORMATION

This MEDICAL DEVICE is intended for use in the electromagnetic environment specified below. The purchaser or user of the MEDICAL DEVICE should ensure that it is used in the specified environment.

Table 1

| Measurement of radio-frequency radiation levels                                                              | Standard    | Characteristics of the electromagnetic environment                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------|-------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Radio interference according to CISPR 11                                                                     | Group 1     | The device uses of radio-frequency energy only for internal functions.<br><br>The level of radio-frequency interference emissions is low and is not likely to disrupt the operation of nearby electronic equipment.                                                                                                                                                                                       |
| Radio interference according to CISPR 11                                                                     | Class A     | The emission characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 Class A). If this equipment is used in residential areas (which typically require CISPR 11 Class B), then adequate operation of the radiocommunication services cannot be guaranteed. The user may need to take mitigation measures such as relocation or reorientation of equipment. |
| Harmonic components of the consumed current according to International Electrotechnical Commission 61000-3-2 | Class A     |                                                                                                                                                                                                                                                                                                                                                                                                           |
| Voltage fluctuations and flicker according to International Electrotechnical Commission 61000-3-3            | Corresponds |                                                                                                                                                                                                                                                                                                                                                                                                           |

This MEDICAL DEVICE is used in a specialized healthcare/radiology setting. The ambient conditions are specified in the INSTRUCTIONS FOR USE.

This MEDICAL DEVICE has been tested for use in a healthcare setting as described above.

However, the high-frequency radiation emission characteristics and immunity may be affected by the connected data cables, depending on the length of these cables and how they are laid.

**This MEDICAL DEVICE is intended for use in the electromagnetic environment specified below. Users of this MEDICAL DEVICE are responsible for compliance with these conditions.**

Table 2

| Guidance and manufacturer's declaration - immunity                                                                                                                                                 |                                                                                                                                                                         |                                                                                                                                                                          |                                                                                                                                                                                           |
|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The MEDICAL DEVICE is intended for use in the electromagnetic environment specified below. The purchaser or user of the MEDICAL DEVICE should ensure that it is used in the specified environment. |                                                                                                                                                                         |                                                                                                                                                                          |                                                                                                                                                                                           |
| Immunity test                                                                                                                                                                                      | International Electrotechnical Commission 60601 test level                                                                                                              | Compliance level                                                                                                                                                         | Electromagnetic environment - guidance                                                                                                                                                    |
| Electrostatic discharge according to International Electrotechnical Commission 61000-4-2                                                                                                           | $\pm 6$ kV - contact discharge<br>$\pm 8$ kV - air discharge                                                                                                            | $\pm 6$ kV - contact discharge<br>$\pm 8$ kV - air discharge                                                                                                             | The floors of the room should be made of wood, concrete or ceramic tiles. If the floors are covered with synthetic material, then the relative humidity of the air should be at least 30% |
| Nanosecond transients according to International Electrotechnical Commission 61000-4-4                                                                                                             | $\pm 2$ kV - for power lines<br>$\pm 1$ kV - for input-output lines                                                                                                     | $\pm 2$ kV - for power lines<br>$\pm 1$ kV - for input-output lines                                                                                                      | The quality of electrical energy in the electrical network of the building must correspond to the typical conditions of healthcare facilities                                             |
| High energy microsecond transients according to International Electrotechnical Commission 61000-4-5                                                                                                | $\pm 1$ kV - when interference is applied according to the "wire-wire" scheme<br>$\pm 2$ kV - when interference is applied according to the "wire-ground" scheme        | $\pm 1$ kV - when interference is applied according to the "wire-wire" scheme<br>$\pm 2$ kV - when interference is applied according to the "wire-ground" scheme         | The quality of electrical energy in the electrical network of the building must correspond to the typical conditions of healthcare facilities                                             |
| Voltage dips, interruptions and voltage variations in input power lines according to International Electrotechnical Commission 61000-4-11                                                          | $<5\%$ $U_t$ (voltage dip $>95\%$ $U_t$ for 0.5 period)<br>$40\%$ $U_t$ (voltage dip $60\%$ $U_t$ ) for five periods<br>$70\%$ $U_t$ (voltage dip $30\%$ $U_t$ ) for 25 | $<<5\%$ $U_t$ (voltage dip $>95\%$ $U_t$ for 0.5 period)<br>$40\%$ $U_t$ (voltage dip $60\%$ $U_t$ ) for five periods<br>$70\%$ $U_t$ (voltage dip $30\%$ $U_t$ ) for 25 | The quality of electrical energy in the electrical network of the building must correspond to the typical conditions of healthcare facilities                                             |

|                                                                                                 |                                                      |                                                      |                                                                                                                         |
|-------------------------------------------------------------------------------------------------|------------------------------------------------------|------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------|
|                                                                                                 | periods<br><br><5% Ut (voltage dip >95%) Ut for 5 s) | periods<br><br><5% Ut (voltage dip >95%) Ut for 5 s) |                                                                                                                         |
| Power frequency magnetic field according to International Electrotechnical Commission 61000-4-8 | 3 A/m                                                | 3 A/m                                                | Power frequency magnetic field levels should be provided in accordance with typical conditions of healthcare facilities |
| Note - Ut - the voltage level of the electrical network before the test action is applied       |                                                      |                                                      |                                                                                                                         |

Table 3

| Guidance and manufacturer's declaration - immunity                                                                                                                                               |                                                                    |                                                                    |                                                                                                                                                                                                                                                                                                                                                                      |
|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------|--------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| The Medical Device is intended for use in the electromagnetic environment specified below. The purchaser or user of the MEDICAL DEVICE must ensure that it is used in the specified environment. |                                                                    |                                                                    |                                                                                                                                                                                                                                                                                                                                                                      |
| Immunity test                                                                                                                                                                                    | International Electrotechnical Commission 60601 test level         | Compliance level                                                   | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                               |
| Conducted disturbances induced by radio frequency electromagnetic fields according to International Electrotechnical Commission 61000-                                                           | 3 V RMS from 150 kHz to 80 MHz<br><br>3 V/m from 80 MHz to 2.5 GHz | 3 V RMS from 150 kHz to 80 MHz<br><br>3 V/m from 80 MHz to 2.5 GHz | The distance between the mobile radiotelephone communication system used and any element of the MEDICAL DEVICE must be at least the recommended separation distance, which is calculated in accordance with the expression below for the frequency of the transmitter. The recommended spacing is:<br>$d = 1,2\sqrt{P}$<br><br>$d = 1,2\sqrt{P}$<br>(80 to 800 MHz); |

4-6

Radiated RF  
electromagnetic field  
according to  
International  
Electrotechnical  
Commission 61000-  
4-3

$$d = 2,3\sqrt{P}$$

(800 MHz to 2.5 GHz).  
where d is the recommended  
spacing, m.

P- rated maximum output power  
of the transmitter, W, set by the  
manufacturer

The field strength during the  
propagation of radio waves from  
stationary radio transmitters,  
based on the results of  
observations of the  
electromagnetic environment,  
should be lower than the  
compliance level in each  
frequency band.

Interference may occur in the  
vicinity of equipment marked  
with the

The field strength of radio waves from fixed radio transmitters, such as radio (cellular/cordless) base stations and land mobile radios, amateur radio, AM and FM broadcast transmitters, television transmitters, cannot be calculated with sufficient accuracy. To do this, practical measurements of the field strength must be made. If the measured values at the location of the MEDICAL DEVICE are greater than the applicable compliance levels, the MEDICAL DEVICE should be monitored to verify that it is functioning properly. If a deviation from normal functioning is detected during monitoring, additional measures must be taken, such as reorientation or relocation of the MEDICAL DEVICE.

Outside the frequency band 150 kHz to 80 MHz, field strengths should be kept below 3 V/m.

#### Notes

1 At 80 and 800 MHz, the higher field strength applies.

2 The above expressions are not applicable in all cases. The propagation of electromagnetic waves is

affected by absorption or reflection from structures, objects and people.

Table 4

| Recommended separation distances between portable and mobile RF communications equipment and the MEDICAL DEVICE                                                                                                                                                                                                                                                                                                                                                                                                                                            |                                                                      |                                                     |                                                          |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------|-----------------------------------------------------|----------------------------------------------------------|
| The MEDICAL DEVICE is intended for use in an electromagnetic environment in which radiated interference levels are monitored. The customer or user of the MEDICAL DEVICE can help prevent electromagnetic interference by maintaining a minimum separation between portable and mobile RF communications equipment (transmitters) and the MEDICAL DEVICE as recommended below, taking into account the maximum output power of the communications equipment.                                                                                               |                                                                      |                                                     |                                                          |
| Rated maximum output power of the transmitter, W                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           | Spatial separation, m, depending on the frequency of the transmitter |                                                     |                                                          |
|                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            | $d = 1,2\sqrt{P}$<br>in the band from 150 kHz to 80 MHz              | $d = 1,2\sqrt{P}$<br>in the band from 80 to 800 MHz | $d = 2,3\sqrt{P}$<br>in the band from 800 MHz to 2.5 GHz |
| 0,01                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       | 0,12                                                                 | 0,12                                                | 0,23                                                     |
| 0,1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 0,38                                                                 | 0,38                                                | 0,73                                                     |
| 1                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          | 1,2                                                                  | 1,2                                                 | 2,3                                                      |
| 10                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         | 3,8                                                                  | 3,8                                                 | 7,3                                                      |
| 100                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        | 12                                                                   | 12                                                  | 23                                                       |
| <p>Notes</p> <p>1 At 80 and 800 MHz, the higher field strength applies.</p> <p>2 The above expressions are not applicable in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.</p> <p>3 When determining the recommended separation distances for transmitters with a rated maximum output power not listed in the table, the expressions given are substituted with the rated maximum output power in watts specified in the transmitter manufacturer's documentation.</p> |                                                                      |                                                     |                                                          |

## RF IMMUNITY OF WIRELESS COMMUNICATION EQUIPMENT

| Frequency range for industrial, scientific and medical equipment (ISM) (MHz) | Operational                                                       | Distance (m) | Immunity test level (W/m) |
|------------------------------------------------------------------------------|-------------------------------------------------------------------|--------------|---------------------------|
| 300-390                                                                      | TETRA 400                                                         | 0,3          | 27                        |
| 430-470                                                                      | GMRS 460; FRS 460                                                 | 0,3          | 28                        |
| 704-787                                                                      | Range LTE 13, 17                                                  | 0,3          | 9                         |
| 800-960                                                                      | GSM 800/900; TETRA 800, IDEN 820; COMA 850; LTE, Range 5          | 0,3          | 28                        |
| 1700-1990                                                                    | GSM 1800; COMA 1900; GSM 1900; DECT; LTE, Range 1, 3, 4, 25; UMTS | 0,3          | 28                        |
| 2400-2570                                                                    | Bluetooth; WLAN; 802.11 b/g/n; RFID 2450; LTE, Range 7            | 0,3          | 28                        |
| 5100-5800                                                                    | WLAN 802.11 a/n                                                   | 0,3          | 9                         |

## EMC (ELECTROMAGNETIC CAPABILITY) PRECAUTIONS

|                                                                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | <p><b>CAUTION:</b></p> <p>Avoid using this equipment when it is installed close to or stacked with other equipment as it may cause malfunctions. If use in such an environment is necessary, this and other equipment should be observed to ensure that it functions properly.</p>                                                                                                                                                                                         |
|    | <p><b>ATTENTION:</b></p> <p>The use of accessories, probes and cables other than those specified or supplied by the manufacturer of this equipment may result in increased levels of electromagnetic radiation or reduced immunity to electromagnetic interference of this equipment, resulting in malfunctions.</p>                                                                                                                                                       |
|    | <p><b>ATTENTION:</b></p> <p>Portable radio communications equipment (including peripherals such as antenna cables and external antennas) should be used at least 30 cm (12 inches) away from any system components, including manufacturer-specified cables. Failure to do so may result in malfunctions of this equipment.</p>                                                                                                                                            |
|   | <p><b>ATTENTION:</b></p> <p>The use of mobile RF communications equipment may affect the operation of the MEDICAL DEVICE. The use of portable or mobile RF devices in the vicinity of any part of the system (including cables) beyond the recommended safety distance is prohibited.</p> <p>The operation of the MEDICAL DEVICE may be affected by communications equipment such as mobile phones. You need to turn it off or put it in airplane mode.</p>                |
|  | <p><b>ATTENTION:</b></p> <p>Cables, converters and other accessories, the use of which as part of the MEDICAL DEVICE ensures compliance with the requirements of GOST R International Electrotechnical Commission 60601-1-2-2014, are an integral part of the system, necessary to ensure its normal operation, cannot be disconnected by the user and were taken into account when measuring EMC (electromagnetic compatibility).</p>                                     |
|  | <p><b>ATTENTION:</b></p> <p>Use of cables, transducers, and accessories not supplied with the MEDICAL DEVICE, with the exception of transducers, accessories, and cables supplied by the system manufacturer as replacement parts for internal components, may result in increased electromagnetic emissions or reduced system immunity</p>                                                                                                                                |
|  | <p><b>ATTENTION:</b></p> <p>The MEDICAL DEVICE should not be used in close proximity or in conjunction with other equipment. When using the product in close proximity to other devices (including when they are placed on top of each other), correct operation is not guaranteed! If the Device for printing is to be used in the vicinity of other devices, a function test should be performed to verify that the system functions properly in this configuration.</p> |
|  | <p><b>ATTENTION:</b></p> <p>All cables used to connect peripherals must be fully shielded and grounded. Use of unshielded and ungrounded cables may result in radio frequency interference.</p> <p>Do not use accessories, transducers or cables other than those supplied with the</p>                                                                                                                                                                                    |

|                                                                                   |                                                                                                                                                                                                                                                                                     |
|-----------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                   | <p><b>MEDICAL DEVICE.</b> An increase in electromagnetic emission or a change in resistance is possible, which can lead to malfunctions of the system! Use only system manufacturer-supplied accessories, transducers, and cables as replacement parts for internal components.</p> |
|  | <p><b>ATTENTION:</b><br/>Do not allow portable radio devices to be used near the MEDICAL DEVICE. Mobile phones near the machine must be turned off. These rules must be observed when the MEDICAL DEVICE is turned on.</p>                                                          |

## CABLES AND ACCESSORIES

Cables and accessories have been tested and found to comply with the related IEC60601-1-2 EMC (ELECTROMAGNETIC COMPATIBILITY) standard:

Attention: The use of cables and accessories other than those specified in this INSTRUCTIONS FOR USE, or other replacement parts than those ordered from the manufacturer, may result in increased levels of electromagnetic interference and/or increased susceptibility to it.

| Function           | Type;<br>Max length. | Comments |
|--------------------|----------------------|----------|
| Network connection | CAT5<br>5 m.         | Shielded |

## APPENDIX A (CODES AND TYPES OF ERRORS)

### Types of Device for printing errors

|                          |                                                                                                                       |
|--------------------------|-----------------------------------------------------------------------------------------------------------------------|
| Back cover error:        | back cover is open                                                                                                    |
| Communication error:     | host software not included                                                                                            |
| Film Press Motor Error:  | The film presser motor cable is not inserted correctly or the main control panel cable is disconnected.               |
| Transmission motor error | the transmission motor cable is not inserted correctly                                                                |
| Print head motor error   | The print head motor cable is not inserted correctly, or the main control panel cable is disconnected.                |
| Display screen error:    | The display screen cable is loose or the main control panel cable is loose.                                           |
| FPGA status error:       | the computer crashes or the USB cable between the motherboard and the computer is in poor contact or is not connected |
| Host error:              | when the device is turned on, the computer host is not turned on, or the application is not running                   |
| Print sensor error       | Check print head motor sensor cable                                                                                   |
| Down Motor Sensor Error: | Check Down Motor Sensor                                                                                               |
| Film size error          | The film does not match the type<br>The film used is not allowed                                                      |

### Codes of errors

| Type                                  | Code       |
|---------------------------------------|------------|
| Tray 1 is opened                      | 0x00000001 |
| Tray 2 is opened                      | 0x00000002 |
| Printing head plug-in error           | 0x00000010 |
| FPGA printing error                   | 0x00000020 |
| Film blocking                         | 0x00000040 |
| Printing head temperature is too high | 0x00000080 |
| Back cover is opened                  | 0x00000100 |
| Communication error                   | 0x00000200 |
| Film Press motor error                | 0x00000400 |
| Motor start error                     | 0x00000800 |
| Print head motor error                | 0x00001000 |
| Display screen error                  | 0x00002000 |
| FPGA status error                     | 0x00004000 |

|                                  |            |
|----------------------------------|------------|
| Host error                       | 0x00008000 |
| Print sensor error               | 0x00010000 |
| Key error                        | 0x00020000 |
| Tray 1 is out of film            | 0x00040000 |
| Film tray error                  | 0x00100000 |
| Tray 2 is out of film            | 0x00080000 |
| Print head Motor Sensor Error    | 0x00200000 |
| Film Pressure Motor Sensor Error | 0x00400000 |
| Film format error                | 0x00800000 |
| No film                          | 0x04000000 |
| RFID chip error                  | 0x08000000 |
| Send sensor error                | 0x10000000 |
| Receive failure                  | 0x20000000 |

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