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CERTIFICATE

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User's Manual

Automatic Chemiluminescence Immunoassay Analyzer
MAGLUMI®, with accessories, in design variants

Dear User! Thank you for choosing automatic chemiluminescence immunoassay analyzer MAGLUMI®.

Please read the manual carefully before using the analyzer to ensure safe and efficient operation of the analyzer.

After reading it, please keep this manual within the convenient reach so that you can refer to it when the need arises.

Should you have any questions related to automatic chemiluminescence immunoassay analyzer MAGLUMI® please contact local company representative office.

Shenzhen New Industries Biomedical Engineering Co., Ltd.
Address: No. 23, Jinxiu East Highway, Pingshan District, 518122
Shenzhen, People's Republic of China
Tel.: +86-755-21536601 Fax: +86-755-28292740

Authorized representative in the RF:

CLS LLC, 125363, city of Moscow,
6 Novoposelkovaya St., bld. 216.1, room
Tel.: +7-495-9469597

Shanghai International Holding Corp. GmbH (Europe)
Address: Eifestrasse 80, 20537 Hamburg, Germany
Tel.: +49-40-2513175 Fax: +49-40-255726

Model	Catalogue number
MAGLUMI 800	23020003
MAGLUMI 2000 Plus	23020007

Intellectual Property Statement

Shenzhen New Industries Biomedical Engineering Co. Ltd (SNIBE) is the intellectual property right holder for this product and the copyright to this manual.

All rights reserved. Unauthorized copying, editing or translation of this user's manual or parts thereof is prohibited.

Snibe, **MAGLUMI**, **Blossays**, **Biolumi**, **Preaccu**, **SnibeLinker**, **SnibeLis**, **SATLARS**, **Molecision**, and etc. are registered trade marks or trade marks owned by SNIBE Co., Ltd. In China and other countries.

Product information

Product name: Automatic chemiluminescence immunoassay analyzer MAGLUMI®, with accessories, in design variants:

Model: Automatic chemiluminescence immunoassay analyzer MAGLUMI 800 with accessories.
Automatic chemiluminescence immunoassay analyzer MAGLUMI 2000 Plus with accessories.

Product: main structural elements are dosing, reaction tracking and data processing systems.

Designation: Automatic chemiluminescence immunoassay analyzer MAGLUMI® with accessories in design variants is used with adapter reagents for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in serum, plasma, urine, whole blood.

User's manual information

Publishing date: 2023-03

Version: 1.1.

Relevant software: 2.14.18.1 and/or compatible version

Table of Contents

INTELLECTUAL PROPERTY STATEMENT	2
PRODUCT INFORMATION	3
USER'S MANUAL INFORMATION	3
TABLE OF CONTENTS.....	1
NOTES.....	1
DESIGNATION.....	2
OPERATING NOTES.....	5
WARNING SIGNS.....	8
OTHER SYMBOLS.....	11
1. ABOUT THIS MANUAL.....	1
1.1 TEXT CONVENTIONS.....	1
1.2 BUTTON	1
1.3 WARNINGS AND NOTES	1
1.4 GLOSSARY.....	2
2. MEASURING PRINCIPLE	1
2.2 MEASURING PRINCIPLE.....	2
2.3 CALIBRATION.....	3
3. SYSTEM DESCRIPTION.....	1
3.1 SYSTEM DESIGN.....	1
3.2 PRODUCT SPECIFICATIONS.....	2
3.3 ANALYZER	3
3.3.3. System liquid.....	12
3.3.4. Waste Liquid.....	12
3.4 COMPUTER SYSTEM.....	22
3.4.1. Computer system components.....	22
3.4.2. Software components	22
3.5 SUPPLIED MATERIALS.....	26
REQUIRED MATERIALS (BUT NOT PROVIDED):.....	26
4. INSTALLATION AND STARTUP	1
4.1 ANALYZER TRANSPORTATION AND STORAGE REQUIREMENTS	1
4.1.1. Transportation and Storage Requirements	1
4.2 INSTALLATION REQUIREMENTS.....	1
4.2.1. Installation Environment Requirements.....	1
4.2.2. Power requirements	1
4.2.3. Space Requirements	2
4.2.4. Temperature and Humidity Requirements	2
4.2.5. Analyzer Transportation and Fastening	2
4.2.6. Language	3
5. DAILY OPERATING PROCESS	1
5.1 DAILY OPERATING PROCESS	1
5.2 TEST PREPARATION	2
5.2.1 Check before Starting.....	2
5.2.2 Power-on and Login the Software.....	2

5.2.3	Check Consumables	3
5.2.4	Analyzer Status Check.....	3
5.2.5	Confirm Test Condition.....	5
5.2.6	Prepare Reagent	5
5.3	TEST ANALYSIS.....	6
5.3.1	Assay Calibration.....	6
5.3.2	Sample Registration.....	9
5.3.3	Start Test.....	10
5.3.4	Additional Samples and Assays.....	10
5.4	TEST RESULTS.....	11
5.4.1	Sample results	11
5.4.2	Results and QC chart.....	11
5.5	END OF ANALYSIS	12
5.5.1	Shutdown.....	12
5.5.2	Operation after System Shutdown.....	12
5.6	[SYSTEM] MENU DESCRIPTION.....	13
5.6.1	<Mode>	13
5.7	[SETTINGS] MENU	15
5.8	[PROCESS] MENU.....	15
5.9	[SYSTEM TEST] MENU.....	16
5.9.1	Load Light Check.....	18
5.9.2	Add Light Check Liquid	18
5.10	[REAGENTS] INTERFACE	19
5.10.1	Reagent Area.....	20
5.11	[REAGENT DATA] MENU	23
5.11.1	Rack station.....	25
5.12	[REPORT] MENU.....	26
6	STATUS DISPLAY AND HANDLING	1
6.1	STATUS DISPLAY	1
6.1.1	Reservoir.....	1
6.1.2	Temperature and Voltage.....	3
6.1.3	Waste.....	3
6.2	RESERVOIR HANDLING.....	5
6.2.1	Placing Cuvette.....	5
6.2.2	Adding Wash Liquid.....	5
6.2.3	Change Starter Reagents	6
6.2.4	Use of cuvette.....	6
6.2.5	Liquid waste.....	7
7	REAGENT LOADING	1
7.1	REAGENT CARTRIDGE.....	1
7.2	REAGENT LOADING	1
7.2.1	Reagent preparation	1
7.2.2	Reagent loading	1
7.2.3	Reagent removal	2
7.3	PROPER REAGENTS STORAGE.....	2
8	SAMPLE LOADING	1
8.1	SAMPLE RACK	1
8.1.1	Sample rack assembly.....	1
8.1.2	Sample rack type.....	2
8.2	LOADING SAMPLE	2
8.2.1	Sample preparation.....	2
8.2.2	Loading sample rack.....	3
8.3	MAINTENANCE OF TEST TUBE RACKS	3

Table of contents

9 SYSTEM MAINTENANCE	1
9.1 ANALYZER MAINTENANCE	1
9.1.1 <i>Daily Maintenance</i>	1
9.1.2 <i>Weekly Maintenance</i>	1
9.1.3 <i>Monthly Maintenance</i>	1
9.1.4 <i>Lasting System Shutdown</i>	2
9.1.5 <i>List of Components</i>	2
9.2 SOFTWARE MAINTENANCE	4
9.2.1 <i>Software Safety Classification</i>	4
9.2.2 <i>Software Update Method</i>	4
10 TROUBLESHOOTING AND DIAGNOSTICS	1
10.1 COMMON FAILURES	1
11 ENVIRONMENTAL PROTECTION REQUIREMENTS	1
12 LIST OF INTERNATIONAL STANDARDS THAT APPLY TO THE MEDICAL PRODUCT	1
13 STORAGE, TRANSPORTATION AND PACKAGING	1
14 MANUFACTURER'S WARRANTY AND LIABILITY	1
APPENDIX A PERFORMANCE INDICATORS	1

Notes

This section contains important information and instructions for safe use of the product.

Please read the following instructions before the use of the analyzer.

Designation

Automatic chemiluminescence immunoassay analyzer MAGLUMI® with accessories in design variants is used with adapter reagents for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in serum, plasma, urine, whole blood.

User's manual is designated for automatic chemiluminescence immunoassay analyzer MAGLUMI®.

This manual describes the design and operation of the analyzer and provides maintenance and troubleshooting recommendations. Make sure to comply with the user's manual when operating the analyzer.

Indications and contraindications for use

Indications:

The MAGLUMI® analyzer is indicated for complex diagnostics and monitoring of specific diseases* and medical conditions.

Used in clinical practice for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in serum, plasma, urine, whole blood.

* **Testing:** including hormones, tumor antigens, proteins and peptides, liver disease, cardiomyopathy, myocardial disease, immune function, autoantibodies, infectious diseases, vitamin and drug concentrations in blood, other enzymes, D-dimer, coagulation, microalbuminuria, human immunoglobulin E, procalcitonin.

Contraindications:

No contraindications to the use of the medical product for in vitro diagnostics and monitoring have been identified.

Safety Note

To ensure safe operation of the system, read the manual carefully before work. Failure to follow safety precautions when operating the analyzer result in personal injury or equipment damage.

The system manufacturer warrants that it meets the safety requirements for electronic and medical analyzers. Legislation establishes requirements for installation and operation of the system that shall be followed by installers and operators.

WARNING!

- 1) Violation of maintenance recommendations may result in malfunction of the analyzer and endanger the health of personnel.
- 2) To ensure safe and reliable operation of the analyzer, installation and maintenance of the analyzer should be performed only by or authorized by our company or authorized personnel and operating engineers. All supplied system components are subject to mandatory testing by our company or authorized dealers.

1. Prevention of injuries caused by moving parts of the analyzer

Familiarize yourself with safety precautions to prevent injuries that can be caused by moving parts of the analyzer in operation.

WARNING!

- 1) Do not touch moving parts of the system or travel ways while the analyzer is operating. Moving parts include the pipette dispenser, incubator, wash station, washing station loading device, incubator loading device or**

- sample loading device, return device, and ejector.
- 2) Do not put objects in the ways of moving parts while the analyzer is in operation. Failure to do so may result in personal injury or analyzer damage.
 - 3) The caps on the sample pipes prevent the pipette aid from moving. It is necessary to remove the caps from the sample pipes.
 - 4) The analyzer is fitted with a lockable lid. Lower the lid and close the lock while the analyzer is in operation. If it is necessary to open the analyzer lid, deenergize power supply to avoid personal injury or equipment damage.

2. Risk of electric shock

Familiarize yourself with safety precautions to prevent electric shock.

WARNING!

- 1) Unauthorized maintenance personnel shall not access the rear and side covers of the energized analyzer.
- 2) Leakage of reagent or sample inside the analyzer may result in damage and electric shock. In such a case, turn off the power supply immediately and contact our technical service department.
- 3) Disconnect power supply before opening the rear and side covers when replacing system components.
- 4) Improper grounding may result in electric shock and damage to the analyzer.
- 5) Ensure that the input voltage complies with the requirements of the analyzer.
- 6) Do not touch or perform electrostatic charge removal on components with static electricity protection warning labels..

3. Fire hazard prevention

Familiarize yourself with safety precautions to prevent fire when working with organic solutions.

WARNING!

- 1) Do not use organic solutions while testing.
- 2) The analyzer is not designed to be explosion-proof. Any organic solutions should be used with caution to avoid the likelihood of fire or explosion.

4. Prevention of the risk of exposure to laser radiation

Familiarize yourself with safety precautions to prevent burns from the barcode scanner laser.

WARNING!

Direct exposure of the barcode scanner laser to eye retina will result in eye injury. Avoid direct eye contact with the scanner's laser ray.

5. Liquid waste disposal

Familiarize yourself with safety precautions to prevent environment contamination and personal injury when disposing of liquid waste.

**WARNING!**

- 1) Certain liquid wastes shall be disposed of in accordance with environmental regulations and toxicity standards. All departments using analyzers shall dispose of liquid waste in accordance with toxicity regulations and the manufacturer's or authorized dealer's recommendations.
- 2) Waste patients' liquids with infectious diseases should be disposed with the use of designated disposal devices.

6. Prevention of chemical exposure hazards

Familiarize yourself with safety precautions to prevent chemical exposure to reagents and consumables.

**WARNING!**

- 1) Please carefully read the safety data sheet for reagents and consumables to familiarize yourself with the safety instructions and precautions for handling the substances.
- 2) Avoid direct contact of reagents and consumables with hands and clothing. In case of accidental contact, wash hands and clothing immediately with soap and water. In case of accidental contact of reagents or consumables with eyes, immediately flush them with plenty of water and consult an ophthalmologist.

7. Prevention of biochemical hazards

Familiarize yourself with safety precautions to prevent biochemical hazard.

**WARNING!**

- 1) Improper handling of samples may result in contamination. Do not allow samples, mixtures or liquid waste to come into contact with your hands or other parts of your body. Always wear gloves and special protective clothing while working; wear protective goggles if necessary.
- 2) Cuvettes contact potentially infected samples, so used cuvettes should be disposed in waste bags to isolate the potential source of infection.
- 3) Handle reagents with care to avoid direct contact of reagents and consumables with hands and clothing. In case of accidental contact, wash hands and clothing immediately with soap and water. In case of accidental contact with eyes, flush them immediately with plenty of water and consult an ophthalmologist.
- 4) If a small amount of sample or reagent gets on the surface of the analyzer, use a soft cloth and alcohol to clean the device. If a large amount of sample or reagent gets on the surface of the analyzer, stop operation immediately and contact an authorized technician to resolve the issue.
- 5) Before transporting the analyzer for a long distance, the device should be thoroughly disinfected to avoid spreading a potential source of infection.

8. Disposal of analyzer

Please follow the below analyzer disposal instructions.

WARNING!

Certain substances in the analyzer shall be disposed in accordance with the

environmental regulations. Disposal shall comply with the local regulatory requirements.

9. Protect your computer from viruses

Familiarize yourself with safety precautions to protect your computer from viruses.

WARNING!

- 1) Use the data transmission and exchange functions only within the permissible limits to protect your computer from viruses or damage to the system software due to improper operation. Sources of computer viruses include floppy disks, USB drives, network, etc.
 - 2) Do not install unauthorized software or hardware that may affect the normal operation of the computer's system software. Do not start other software while the system is running.
-

Operating Notes

Please carefully read the following operating precautions for proper and efficient operation of the analyzer.

1. General precautions

Use of this analyzer requires an understanding of its operation and general safety precautions applicable thereto. Failure to follow the instructions specified in this manual may adversely affect the protection of the analyzer.

NOTE

- 1) For clinical evaluation, use the results of the assay in conjunction with clinical symptoms or the results of other investigations.
 - 2) Operating instructions are subject to revision without advance notice. If necessary, please contact a customer service representative.
 - 3) This analyzer is designed for professional use by medical laboratory personnel, qualified physicians, nurses and laboratory technicians only.
 - 4) When working with chemicals, do not touch the computer monitor, mouse, or keyboard.
 - 5) Do not bend or squeeze the drain pipe to avoid making it difficult to drain the liquid and avoid its spillage, as in extreme cases this may result in product failure.
 - 6) The analyzer emits a large amount of heat during operation, which is emitted from the backside of the product. Therefore, the room, where the analyzer is installed should be well ventilated to ensure that the analyzer is cooled and ventilated if necessary. Make sure that the airflow is not directed directly at the analyzer, as this may affect the quality of the results.
 - 7) Before the first use, adjust all analyzer components so that they operate accurately.
 - 8) Make sure that there are no bubbles in the starting reagent and washing liquid. Otherwise, the correctness of the analysis results cannot be guaranteed.
 - 9) Check the correct installation of the starting reagent containers; do not use mixed starting reagents. Otherwise, the correctness of the analytical
-

results cannot be guaranteed.

- 10) Use new or cleaned cuvettes to ensure safe operation of the analyzer and accurate assay results.
- 11) To ensure safe analyzer operation and stable assay results, do not use expired wash liquid.
- 12) The analyzer should be turned on at least 30 minutes before use to ensure stabilization of the measurement system. The accuracy of the incubator temperature setpoint should be within $\pm 0.5^{\circ}\text{C}$ of the setpoint, and its deviation should not exceed 0.3°C .
- 13) Prior to analysis ensure that consumables (wash liquid, starting reagents, cuvettes, reference solution, etc.) meet the requirements.
- 14) Perform a background wash (BGW) at least once before analysis to ensure that the results of this test are within the normal range. Otherwise, the correctness of the assay results cannot be guaranteed.
- 15) No bubbles should appear on the sample surface during pipetting, as this may cause a pipetting error. Do not remove the reagent until the analyzer finishes the pipetting process.
- 16) Quality control shall be performed when using this product to analyze samples. Otherwise, the correctness of the assay results cannot be guaranteed.
- 17) The analyzer is fitted with the error indicator. If an error occurs in the analyzer, the indicator lights up and an alarm sounds. In such a case, clear the error to ensure normal analyzer operation and accurate assay results. Please refer to Chapter 11 to learn more about the troubleshooting methods.
- 18) Use the recommended solution to clean the system pipes and keep them in functional condition.
- 19) Any serious adverse event related to medical products shall be reported to the manufacturer as well as to a designated employee.
- 20) If a hazardous substance gets on the surface or inside the product, appropriate disinfection should be applied; do not use cleaners or sanitizers that may chemically react with the components of the products or materials contained inside the product and cause a hazard; if there is any doubt about the compatibility of sanitizers or cleaners with the components of the product or materials contained therein, please contact the manufacturer or its authorized representative in the Russian Federation.

2. Operation environment

NOTE

This analyzer shall be installed in operation environment that meets the requirements of this manual. Failure to install the analyzer under the specified conditions may result in inaccurate results and damage to the system.

3. Electromagnetic compatibility

NOTE

- 1) Automatic chemiluminescence immunoassay analyzer MAGLUMI® (MAGLUMI 2000 Plus, MAGLUMI 800) with the accessories complies with the requirements for immunity and emissions of IEC 61326-2-6.
- 2) It is the responsibility of the user to ensure electromagnetic compatibility

conditions for proper operation of the analyzer.

- 3) It is advised that the electromagnetic conditions be evaluated before the use of the analyzer.
-

WARNING!

- 1) 1)Automatic chemiluminescence immunoassay analyzer MAGLUMI® (MAGLUMI 2000 Plus, MAGLUMI 800) with the accessories has been designed and tested to meet the requirements of IEC/CISPR 11 for Class A analyzers. This analyzer may cause radio interference in residential areas and precautions should be taken to avoid such interference.
 - 2) Do not operate the analyzer near strong radiation sources (e.g., unshielded RF sources) as this may interfere with normal operation of the product.
 - 3) This product is not intended for use in residential accommodation.
-

4. System maintenance

ATTENTION!

- 1) Follow the instructions specified in this manual when performing regular maintenance of the analyzer. Improper maintenance can adversely affect the accuracy and reproducibility of assay results and can result in malfunction or failure of the analyzer.
 - 2) Disconnect all power supplies to the system and unplug the power connector before maintenance and repair. Failure to do so may result in malfunctions or damage the analyzer.
 - 3) Traces of potentially infected samples may remain on the surfaces of the analyzer. Always wear gloves and special protective clothing to protect against infection when performing maintenance and repair.
 - 4) Dust may form on the surface of the analyzer during long-term storage. Use a soft, wet cloth to clean the surface. Do not allow penetration of water drops into the analyzer.
 - 5) The analyzer has no user-serviceable parts. Do not attempt to disassemble the analyzer housing or components. If necessary, contact the authorized company personnel for assistance.
-

5. Samples, reagents and controls

WARNING!

- 1) The presence of drugs, anticoagulants and preservatives in samples may affect the results of some assays.
 - 2) Ensure proper sample storage. Improper storage can lead to a change in sample composition and consequently affect the assay results.
 - 3) To prevent evaporation, do not expose samples to air for a long time. Sample evaporation may affect the assay results.
 - 4) Improper storage of reagents and reference solutions, even if not expired, may result in inaccurate assay results and poor system performance. Follow the manufacturer's instructions when using reagents and reference solutions.
 - 5) Calibration is required right after the reagents are changed. The assay results may be inaccurate without calibration.
-

6. Data backup

WARNING!

Store system data on a hard disk of your PC. However, if data is deleted from the hard disk or if it is damaged due to certain circumstances, data recovery will not be possible. It is necessary to regularly backup the assay results and parameters and save them on external media such as CD-ROMs.

Warning signs

Warning infection (WARNING!)

This sign refers to a biological hazard. It can be found on zones of the analyzer where there is a risk of biological contamination, including:

- the frontage of the waste container;
- the frontage of the waste tank;
- surface of the sample area;
- surface of the reagent area.

No Mixing

This sign is located in the area where the solution is placed to remind not mixing up the solution together. It is on the tray of the starter.

Warning danger (Attention! Hazard)

This sign is located in the area where it is easy to get hurt to remind the safety. It is on below the hinge in the middle of the main support, as well as:

- interior the reagent area
- interior of the sample area.

Watch out for the laser

This sign is located in the area with laser beam to the danger of the laser beam. It is on the shell of the product.

Laser window

This sign is located at the laser beam exit window. It is in the right side of the interior of the sample area

Watch out for the movement of moving component

This sign is located in the moving part of the analyzer to remind not touching the moving component during operation. It is in the frontage of the pipetting arm.

Do not actuate
during operation

Watch out for your safety when opening the cover

This sign is to remind not opening the cover when the analyzer is working. It is on the positive side of the handle on the cover.

Don't open the
cover when analyzer
working

Warning danger

Warning hands pinching

This sign is located on the component with squeezing moving part to remind the danger of clamping hand. It is on the plate covering the incubator stacker pipetting area.

Warning hands pinching

Warning electric shock

Pay attention the words on the sign. It is on the top right side of the shell on the back side of the analyzer.

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Other symbols


Symbol	Description
	Manufacturer
	Manufacturing date
	Catalogue Number
	In Vitro diagnostic medical product
	Serial Number
	Authorized Representative in the European Community
	Attention. Please refer to the documents attached

This symbol means that this product should not be disposed with household waste. Proper disposal of this product will prevent environmental and health impacts that may be caused by hazardous substances present in the product. For more information, please contact the distributor who supplied this product.

This side up

This sign is to remind the direction of the package should be upright during transportation. It is on the frontage of the package
the frontage of the wooden box


Keep away from rain


This sign is to remind keeping the package from rain during transportation. It is on the frontage of the package
the frontage of the wooden box.

**WARNING! Fragile.**

This sign is to remind fragile subject inside, handle with care. It is on the frontage of the package
the frontage of the wooden box

**Rolling is forbidden**

This sign is to remind not rolling the package during transportation. It is on the frontage of the package.
the frontage of the wooden box.

**Main switch****Main Switch****F6AL 250V****Submain switch****Submain
Switch****Energized (on)****Deenergized (off)**

1. About this manual


1.1 Text conventions

To facilitate quick understanding and use of this manual, text styles occur in this document are defined as follows:

- Menu items, interface and dialog box names are boldfaced and put in the symbol []. For example, [Settings] menu, [Test] interface and [Assay user settings] dialog.
- Button names are boldfaced and put in <>. For example, <OK> and <Add> buttons.
- User input is boldfaced and appears in “”. For example, [Sample Pipetting Volume] “2” [μl].

1.2 Button

Common buttons names

Button	Description
	Red means the assay is not selected.
	Green means the assay is selected.
	Next page.
	Previous page.
	Move to next page.
	Page up.
	Move to previous line.
	Move to next line.
	Page down.
	Move to last page.

1.3 Warnings and notes

Sign	Words	Description
	WARNING!	Read the statement next to the symbol. The statement is alerting you to an operating

		hazard that can cause personal injury.
	ATTENTION!	Read the statement next to the symbol. The statement is alerting you to a possibility of system damage or unreliable results. \
	NOTE	Read the statement next to the symbol. The statement is alerting you to information that requires your attention.


1.4 Glossary

Glossary	Description
Analyzer	The device, but not PC, printer and connection cables.
Reverse device	In second pipetting step, transfers cuvettes to proper positions in the incubator.
Barcode Reader	Device to read the sample barcode.
BGW	Background wash, test to check quality of analyzer washing.
Chamber	Measurement unit.
Cuvette	Every cuvette has six cavities plastic module, in which immunometrical reaction can take place.
CV%	Coefficient of variation, shows dispersion rate of measurements.
Incubator	Assembly, where cuvettes are incubated and pipetted.
LC-le	Light check of left pipettor is used to test accuracy of pump volume and stability of the analyzer.
LC-ri	Light check of right pipettor is used to test accuracy of pump volume and stability of the analyzer.
Loader	Includes incubator loader and washer loader.
Pipettor	Device for pipetting batches and samples for assay.
Pump system	For high-precision pipetting, washing and starter injection.
Pusher	Exchanges cuvettes between washer and chamber..
Reader	Includes barcode reader and RFID reader.
Reagent Area	Loading area for reagent.
RFID	Radio Frequency Identification Devices, micro-chip present on reagent cartridges to allow recognition and data storage.
OCE	Relative Light Unit (signal measurement unit).
Samples	Anything that can be introduced by operator into sample area racks, including patient liquid samples, reference solutions and external calibrators.
Samples Area	Loading area for samples.
Samples Rack	12 positions module to host sample pipes for MAGLUMI 2000 Plus. 8 positions module for MAGLUMI 800 analyzers.
Stacker	Unit to store cuvettes.
Starter Reagents	Reagents dispensed during the reading to generate chemiluminescent signal.
Washing liquid	Solution to wash pipetting needles and magnetic microbeads after reaction. It is prepared by dilution of concentrate for cleaning with distilled water.
Washer	For washing unreacted material in cuvettes.
Waste Bag	Container for used cuvettes.


2. Measuring Principle

2.1 Assay Procedure

One step assay

Two step assay

Notes


Figure 4-1. Assay procedure

2.2 Measuring Principle

The analyzer uses the principle of immune response by combining inhomogeneous magnetic separation with chemical luminescence, using an ABEI light-emitting marker directly to mark an antigen or antibody. It reacts with antibodies or antigens in the sample, as well as with antigens or antibodies coated with magnetic microbeads, forming an immune antigen-antibody complex. After removing the free antigen or antibody and other impurities using a magnetic field, a substrate is added to the antigen-antibody immune complex, which lights up the marker. By measuring the luminescence intensity and based on the working curve, the concentration of the analyte is determined.

The analyzer's photomultiplier is used to detect light produced in chemiluminescence reaction, within the wavelengths range between 300nm to 650nm. The light peak of the chemiluminescence is emitted at a wavelength of 420nm. The light produced in hemiluminescence reaction is emitted to the photomultiplier and reaches the photocathode plane through the incidence window, triggering photons on the photocathode plane emitting photoelectrons into a vacuum. Photoelectrons accumulate at the first dynode through the focusing electrode, pass subsequent dynodes for secondary electron multiplication, and then secondary electrons emitted from the last dynode are output through the anode. The photomultiplier anode collects secondary electrons after multiplication by dynodes and outputs current signals through an external circuit.

To eliminate differences between photomultipliers and ensure test result consistency between different analyzer, relative light unit (RLU) is used as the unit of measurement for original data.

After being pipetted to the cuvette, samples and reagents are blended, washed and separated before the cuvette is sent to the chamber. Starter 1 is injected into the first hole in the cuvette bar, and then Starter 2 is injected into the same hole after 2.5 seconds, triggering chemiluminescence reaction. Detection of optical signals starts 0.1 second after the chemiluminescence reaction and obtains optical signals of 3.0 seconds. Repeat this step to detect the other five holes in the cuvette bar.


Figure 2-1. Kinetic curve

2.3 Calibration

Because there are differences between the actual working environment and laboratory environment, the master curve should be adjusted to generate the working curve that meets the actual work environment.

Brief description:

- The master curve is determined by 10 standard calibrators.
- Compare two calibration RLUs obtained by calibrators with RLUs of related concentration on the master curve.
- Calculate the difference between two calibration RLUs obtained by calibrators and RLUs of related concentration on the master curve, and carry out linear inference using the recalculated RLU (Y-axis) and concentration (X-axis).
- Calculate RLU differences of other calibrators on the master curve with the correction curve and recalculate the RLU and concentration.
- The recalculated curve is a valid working curve (see Figure 2.3-1. Calibration principle).


Figure 3-1. Calibration principle

3. System Description

Automatic chemiluminescence immunoassay analyzer MAGLUMI® with accessories in design variants and sets of supporting diagnostic reagents constitute a precise trace assay system, which features direct chemiluminescence immunoassay based on magnetic separation of ABEI markers. It is used to quantitatively or qualitatively analyze analytes in blood serum, blood plasma, urine, whole blood. The analyzer automatically completes sample and reagent pipetting, incubation, washing, measurement and result calculation, reducing errors in assay results and improving their accuracy and reproducibility.

3.1 System design

The analyzer is composed of material supply module, tubing system module, temperature control module, mechanical transmission module, optical path detection module, circuit control module and waste collection module.

- The material supply module includes stacker, sample area module, reagent area module, system liquid module, waste liquid module, starter module.
- The tubing system module consists of a pipetting tubing system, system liquid tubing system, optical path detection tubing system and condensate water tubing system.
- The temperature control module consists of an incubator heating module, incubator loader heating module, back transport pipetting area heating module, sample cooling module, reagent cooling module and photomultiplier cooling module.
- The mechanical transmission module consists of cuvette loader, stacker, reagent shaker, incubator, incubator loader, washer loader, washer transport, washer lift, back transport, pusher, chamber transport, chamber lift and pipettor.
- The optical path detection module consists of a chamber module, photomultiplier module and main control circuit.
- The circuit control module consists of a power module, main control board, wire harness and a variety of sensors and motors. The computer is an optional accessory.
- Waste collection module consists of a waste bag, cuvette waste container, polyethylene canister.

3.2 Product specifications

Item		Description
Main function	Test speed	MAGLUMI 2000 Plus, MAGLUMI 800 : 180 tests/hpur, first result delivery time — 17 min. 24-hours standby mode
	Samples area	MAGLUMI 800: accommodates 40 samples MAGLUMI 2000 Plus accomodates 144 samples 1. Continuous loading and cooling function. 2. Temperature: 8 °C ~ 18 °C
	Reagents area	MAGLUMI 800: 9 channels MAGLUMI 2000 Plus 25 channels 1. Continuous loading and cooling function. 2. RFID mark stores information. 3. Easy and fast adding and updating. 4. Temperature: 5 °C ~ 15 °C
	Testing modes	Random sampling mode, sequential sampling mode, STAT
	Sample types	Blood serum, plasma, whole blood, urine
	Barcode type	Supports five barcode formats: Code 128, Code 39, Code 93, Codabar and 2/5 Interleaved
	Sample identification method	Barcode reading or manual input
	Quality control	Bag QC, monthly QC
	Testing specifications	1. Continuous loading, perform STAT procedure first 2. Automatic dilution or selection of dilution degree 3. Aggregate analysis according to specified parameters 4. Automatic receipt of sample data after connection to LIS system 4. 5. Two-point master curve correction method
	Analysis method	Competition method, sandwich method, with specific grip, indirect
	Quantity of cuvettes to be loaded	MAGLUMI 800: 40 cuvette strips MAGLUMI 2000 Plus 120 cuvette strips
	Sample probe	Liquid level measurement, clot detection, pipe collision prevention, liquid level measurement
	Reagents probe	Liquid level measurement, pipe collision prevention, liquid level measurement
	Incubator	MAGLUMI 800: incubation of 13 cuvette strips simultaneously (78 samples) MAGLUMI 2000 Plus: incubation of 16 cuvette strips simultaneously (96 samples) Temperature: 36.8 °C ± 0.5 °C
	Wash station	1. Efficient wash station with overflow protection. 2. Intelligent magnetic field design ensures good separation
External interface	Chamber	1. High-sensitivity photomultiplier pipe with low noise level 2. Reaction module tracking function; and overflow protection function 3. Single photon counting technology
	Display	Contrast: greater than or equal to 1000: 1 Power: ≤ 60 Wt 100 V - 240 V Grid frequency: 50/60 Hz Interfaces: 1 * DisplayPort, 1 * HDMI
	Printer	Supports various common printer brands

System description

	LIS system	LIS system support for two-way communication via ASTM protocol
	Operation system	Supports Windows
Operating conditions	Temperature	10~30 °C
	Relative humidity	≤70 %
	Atmospheric pressure	85.0-106.0 kPa
	Other	Do not place the device near of electromagnetic field sources that may interfere with its operation
Storage conditions	Temperature	-20 ~ 55 °C
	Relative humidity	≤93 %
	Atmospheric pressure	50.0-106.0 kPa
	Other	The room shall not allow penetration of direct sunlight and corrosive gases; the room shall be well ventilated
Safety classification	Safety class	Class I (PE connected)
	Overvoltage reliability category	II
	Degree of contamination	2
Analyzer main unit	Dimensions Length (left to right) x width (front to back) x height, weight	MAGLUMI 800: 102 x 72 x 56 cm 73 kg MAGLUMI 2000 Plus: 144 x 76 x 152 cm 228 kg
	External package dimensions Length (left to right) x width (front to back) x height	MAGLUMI 800: 102 x 80 x 78.4 cm MAGLUMI 2000 Plus: 176 x 94 x 178.2 cm
	Power source	100–240 V AC, 50/60 Hz
	Consumed power (VA)	MAGLUMI 800: 630 B·A MAGLUMI 2000 Plus: 840 B·A
	other	Continuous working hours: 24 hours MTBF: one year
	Dimensional and weight tolerance is ±5% unless otherwise specified.	

3.3 Analyzer

Design (basic configuration) and accessory list (full range of medical product options) and technical specifications

Material supply module


Figure 3-1. Description of material supply module of MAGLUMI 800


Figure 3-2. Components of material supply module of MAGLUMI 2000 Plus

- | | | |
|------------------------|-----------------------|-------------------------|
| 1. Stacker | 2. Sample area module | 3. Barcode reader |
| 4. Reagent area module | 5. Starter module | 6. System liquid module |

Pipe system module


Figure 3-1. Components of pipe system module of MAGLUMI 800


Figure 3-2. Components of pipe system module of MAGLUMI 2000 Plus

- | | |
|-------------------------------------|----------------------------------|
| 1. Dosing pipe system module | 2. Liquid pipe system |
| 3. Optical path pipe system | 4. Condensate pipe system |

Temperature control module


Figure 3-3. Components of temperature control module of MAGLUMI 800


Figure 3-4. Components of temperature control module of MAGLUMI 2000 Plus

- | | |
|---|-----------------------------------|
| 1. Incubator loader heating module | 2. Incubator heating module |
| 3. Back transport pipetting area heating module | 4. Sample cooling module |
| 5. Reagent cooling module | 6. Photomultiplier cooling module |

Mechanical transmission module


Figure 3-5. Components of mechanical transmission module of MAGLUMI 800


Figure 3-6. Components of mechanical transmission module of MAGLUMI 2000 Plus

- | | | |
|--------------------|-----------------------|----------------------|
| 1. Cuvette loader | 2. Stacker | 3. Incubator loader |
| 4. Washer loader | 5. Incubator | 6. Washer transport, |
| 7. Washer lift | 8. Dispenser | 9. Back transport |
| 10. Pusher | 11. Chamber transport | 12. Chamber lift |
| 13. Reagent shaker | | |

Optical path detection module


Figure 3-7. Components of optical path detection module of MAGLUMI 800


Figure 3-8. Components of optical path detection module of MAGLUMI 2000 Plus

- | | | |
|-------------------------|-------------------|---------------------------|
| 1. Main control circuit | 2. Chamber module | 3. Photomultiplier module |
|-------------------------|-------------------|---------------------------|

Circuit control module


Figure 3-9. Components of circuit control module of MAGLUMI 800


Figure 3-10. Components of circuit control module of MAGLUMI 2000 Plus

- | | |
|-------------------------------|-------------------------|
| 1. Power module | 2. Wiring harness |
| 3. Various sensors and motors | 4. Main control circuit |

*NOTE: PC is an optional accessory.

Waste collection module


Figure 3-11. Components of waste collection module of MAGLUMI 800


Figure 3-12. Components of waste collection module of MAGLUMI 2000 Plus

1. Polyethylene canister 2. Cuvette waste container 3. Waste bag

3.3.1. Barcode reader

WARNING!

A barcode reader emits laser beams that are harmful to eyes. Therefore, do not look into the barcode reader.

The barcode reader is located between the sample area and reagent area. When you open the sample area door, the barcode reader starts up automatically.

After a sample rack is inserted, barcode labels of the rack and sample pipes are automatically read. The inserted rack is displayed in [Sample area] interface, and sample information is displayed in **Sample Info**.

Table 3-15. Printing Requirements for Sample Pipe Barcodes

Supported Coding type	Data Length Range (characters)	Parity	Barcode Width (mm)	Barcode Height	Recommended Width (mm)	Recommended Height (mm)
Code128	1–25	Mandatory	0.3–0.8	N/A	0.33	10
Code39	1–25	Optional	0.3–0.8	N/A	0.33	10
Codabar	1–2	Optional	0.3–0.8	N/A	0.33	10
Code93	1–25	Mandatory	0.3–0.8	N/A	0.33	10
2/5 Interleaved	2–24	Optional	0.3–0.8	N/A	0.33	10

Blank area at both sides shall be at least 7 times the width of the barcode.

3.3.2. Starter

Starter reagents containers are located on the right of the analyzer. S1 and S2 are marked to indicate Starter 1 and Starter 2, respectively. Starter status is detected by a liquid level detector in the starter container. After starter supplement or replacement, you need to click **<SystemTest>** in the menu bar to prime the tubing system to ensure that the tubing system is filled with starter.

CAUTION!

Do not spill starter in this area!

3.3.3. System liquid

System liquid is used to wash pipetting needles and magnetic microbeads after reaction.

The system liquid container has a liquid level detector for detecting remaining liquid amount. After supplement or replacement of system liquid, you need to click <SystemTest> in the menu bar to fill up the tubing system.

The ports are connected to the analyzer using a coupling head. Press the metal clip to release and remove the pipe.

3.3.4. Waste Liquid

Sources of liquid waste include magnetic microparticles, starters, system liquid, wash station, and cuvette liquid (patient samples and laboratory reagents).

WARNING!

Biological waste must be disposed of in accordance with related laboratory regulations.
Do wear protective gloves!



Waste bag and waste bin status is displayed by <Waste Status> icon.



3.3.5. Fuses


Location of fuses 3A, 5A, 6A, 10A, 15A.


For MAGLUMI 800:


Fuse specification	Installation	Component arrangement	Number
3A	F2 fuse holder on PCB M800-02-E00	Analyzer Rear Bracket	1, 8
5A	F1 fuse holder on PCB M800-05-E02	Rear side of the analyzer reagent area	2, 9
	F1 fuse holder on PCB M800-05-E00	Analyzer Rear Bracket	3, 8
	F1 fuse holder on PCB M800-02-E00	Analyzer Rear Bracket	1, 8
	F1 fuse holder on PCB M800-03-E00	Analyzer Rear Bracket	4, 8
	F1 fuse holder on PCB M800-01-E00	Analyzer Rear Bracket	5, 8
	F1 fuse holder on PCB M600-04-E00	Analyzer Rear Bracket	6, 8
6A	Filter	Analyzer main switch	7
10A	Not used	-	-
15A	Not used	-	-











For MAGLUMI 2000 Plus:


Fuse specification	Installation	Component arrangement	Number
3A	F1 fuse holder on PCB M2000-05-E00	Analyzer main control board bracket	1
	F1, F2 fuse holder on PCB M2000-06-E00	Analyzer main control board bracket	2
	F1, F2 fuse holder on PCB M2000-08-E00	Analyzer main control board bracket	3
5A	F1 fuse holder on PCB M2000P-05-E01	Rear side of the sample area	11,12
	F1 fuse holder on PCB M2000-09-E00	Analyzer main control board bracket	4
	F1 fuse holder on PCB 06-E01B (M2000P-06-E01)	Rear side of the reagent area	5,10
	F1 fuse holder on PCB M2000-07-E01	Top side of the analyzer main control board bracket	6,9
6A	Filter	Main switch on the left side of the analyzer	7
10A	F3 fuse holder on PCB M2000-00-E00	Front side of the analyzer main control board bracket	8,9
15A	F1, F2 fuse holder on PCB M2000-00-E00		

3.4 Computer system

A computer system is a PC with a preinstalled operating analyzer software. It is composed of a monitor, computer, keyboard and mouse.

- Minimum computer configuration: CPU frequency ≥ 2.0 GHz, hard disk ≥ 240 GB, memory ≥ 4 Gb; RS-232 interface; USB interface; RJ45 network interface; monitor; keyboard and mouse.
- Supported operating systems: Windows.
- Software Data Interface: Connect to the analyzer via a wired network interface.
- User Access: provides user management functions. It supports rights management and is designed for two types of users: administrators and users.
- Network Conditions: The software is connected to the analyzer via a wired network (TCP/IP, RS-232), and the analyzer does not need to be connected to a local area network (LAN).

3.4.1. Computer system components

Operating software is installed on the computer to control operation and data processing. The computer system consists of a computer, LCD monitor, keyboard and mouse.

3.4.2. Software components


Figure 4-1. Software interface

The software has 5 different levels of windows:

1. Menu Bar: displays different menu. It has 9 menu buttons.


figure 4-2. Menu bar

Button	Description
<Start>	It is for starting measuring samples or controls.
<System test>	It is for entering system test functions e.g. flush pipe system and carry out internal test measurements.
<Sample area>	Switch to the [Sample area] interface to perform sample registration.
<Reagents>	Switch to the [Reagents] interface to perform calibration and etc.
<Process>	It is for entering process functions e.g. automatic clearance of the cuvettes.
<Settings>	Customizing assays, controls, dilutions, groups and classes of assays
<System>	Displays information about the system, system settings, maintenance, etc.
<Report>	It is for managing results.
<Home>	Return to home.

2. Interface: displays function interface in different menu, including [Home] interface, [Sample area] interface, [Reagents] interface and submenu interface.

Figure 4-3. [Home] interface

3. Status Bar: display analyzer status and relogin or exit the software.


Figure 4-4. Status bar

Icon	Button	Description
	< Message Box>	It is used for showing the actual system and error messages.

	< Pause Measure>	It is for emergency stop of the analyzer..
	< Reservoir Status>	It indicates the status of the system liquid, starter reagents and cuvettes.
	<Temperature/voltage>	It indicates the supply voltage and temperature.
	<Waste status>	It indicates the status of waste liquid and waste cuvette.
	< Relogin>	It is used for relogin the software.
	< Close Software>	It is used for exiting the software.

The following information is displayed at the bottom of the status bar:

- Currently authorized users and logon status. Possible logon statuses are: "Logon completed", "Logout completed", "Invalid password", "Invalid user name".
- COM Status Indicator. Green indicator indicates that the COM system has been successfully opened, red indicates that COM cannot be opened.
- PLC status indicator. Green LED indicates successful PLC connection, red LED indicates connection to PLC has been failed.

- **4. Submenu:** click the **<Process>**, **<Settings>**, **<System>** or **<Report>** button on the menu bar, on the left display its submenu buttons.


Figure 4-5. Submenu [System]

5. Dialog

A dialog box can open one or more sub-dialog boxes.

There are two different types of dialog boxes.

- a) Clicking submenu button triggers a dialog. For example, click **<Initialize>** button in **[Process]** menu to open message confirmation dialog.


Figure 4-6. Message Dialog

- a) Clicking a button in submenu interface triggers a dialog. For example, click <Edit> button in [Test] interface to open [User Specific Assay Data] dialog.

Figure 4-7. Dialog [User Specific Assay Data]

3.5 Supplied materials

See packing list.

Required materials (but not provided):

Material	Used models	Catalog numbers	SWG No.
In vitro reagent kit for automatic analyzers MAGLUMI®	All modes	individual	individual
(Reaction Modules for automatic analyzers MAGLUMI®	All models	630003	P3H 2021/14862
Starter kit (Starter 1, Starter 2) for automatic analyzers MAGLUMI®	All models	130299004M	P3H 2021/14858
Light Check for automatic analyzers MAGLUMI®	All models	130299006M	P3H 2021/14861
Wash Concentrate for automatic analyzers MAGLUMI®	All models	130299005M	P3H 2021/14863
System Tubing Cleaning Solution for automatic analyzers MAGLUMI®	All models	130299007M	P3H 2021/14881
Waste Bag for automatic analyzers MAGLUMI®	All models	21060726	P3H 2021/14988

4. Installation and Startup

This chapter describes how to install and start Maglumi series analyzer.
Basic analyzer installation is implemented by engineers trained and authorized by Snibe.

NOTE

- 1) To ensure user safety, installation and startup of analyzer can only be implemented by professionals trained and authorized by SNIBE Co., Ltd.
- 2) Analyzer must run under the operating conditions specified in this manual.

4.1 Analyzer Transportation and Storage Requirements

4.1.1. Transportation and Storage Requirements

- The analyzer must be transported in the upright direction.
- Analyzer must be free of moisture, water, violent vibration and extrusion during transportation, and handled with care during loading and unloading.
- The analyzer should be stored in a well-ventilated indoor environment without chemicals, corrosive gases or strong sunlight.
- Store and transport at the temperature $-20^{\circ}\text{C} \sim 55^{\circ}\text{C}$ and relative humidity $\leq 93\%$.
- Atmospheric pressure during transportation and storage 500 hPa $\sim$ 1060 hPa

Disconnect the power cable from the circuit before storage. Protect the analyzer from dust and other contaminants during storage. It is recommended to store the analyzer in its original package. The analyzer should be individually packed before transportation.

4.2 Installation Requirements

4.2.1. Installation Environment Requirements

- For indoor installation and use only. The installation environment should be well ventilated and free of dust, mechanical vibration, loud noise and power interference.
- The table should be flat and able to withstand a minimum weight of 73 kg while installing MAGLUMI 800 ($\pm 5\%$, unless otherwise specified).
- The floor should be flat and able to withstand a minimum weight of 228 kg while installing MAGLUMI 2000 ($\pm 5\%$, unless otherwise specified).
- When the analyzer is working properly, the highest volume at 1 m distance is 40 db.
- No corrosive or flammable gases.

4.2.2. Power requirements

- Voltage: 100-240 V AC, frequency: 50/60 Hz
- Power

Model	Nominal power	Applicable type of power supply cable
MAGLUMI 800	630 Вт	250 V AC, 10 A, cable length ≤ 3 m

MAGLUMI 2000 Plus	840 W	250 V AC, 10 A, cable length $\leq$ 3 m
-------------------	-------	--

- Circuit breaker

Model	Circuit breaker type
MAGLUMI 800	250 V, 6 A, quick grounding
MAGLUMI 2000 Plus	250 V, 6 A, quick grounding

- To ensure reliable operation of the analyzer, it shall be connected to a properly grounded outlet. The outlet shall be conveniently located.

WARNING!

Failure to comply with the above requirements may affect the performance of the analyzer, cause damage to it, or even result in personal injury.

4.2.3.Space Requirements

To ensure space required for repair and maintenance by the personnel, analyzer installation must meet the following requirements:

- The distance between the analyzer rear and the wall cannot be less than 50 cm.
- The distance between left and right sides of the analyzer and the wall cannot be less than 50 cm.
- The distance between the analyzer front and other analyzer cannot be less than 100 cm
- The power outlet should be located in a place with enough space to facilitate connecting and disconnecting the power cable. Do not place the analyzer in a place where it is difficult to operate the disconnecting device.

4.2.4.Temperature and Humidity Requirements

- Ambient temperature: 10°C~30°C, air-conditioner recommended;
- Relative humidity: $\leq$ 70% %

4.2.5.Analyzer Transportation and Fastening

- The analyzer should be kept in an upright position all the time during transportation.
- Minimize vibration during transportation. After transportation, check and commission the analyzer before use.
- Adjust foot height to ensure analyzer levelness when fastening the analyzer.

4.2.6. Language

The standard software is in English. You can click the **<System>** button on the main menu to open the **[Settings]** interface and set the desired language on the “Language” tab.


Figure 4.2-1 [Settings] interface

4.2.7. Performing System Test

When analyzer initialization is complete and there are no error messages or pop-up windows, perform system testing as follows:

- 1) Perform the component prime test according to the below requirements.

Table 2-1. System Test (round 1)

		Number of Times
Priming	Pipettor	3
	Washer	6
	Chambers	3
Cuvettes	BGW	0
	LC-le	0
	LC-ri	0

- 1) After the prime test is complete, perform the background test and left and right pipetting needle checks according to the requirements.

Table 4.2 2. System Test (round 2)Таблица

		Number of Times
Priming	Pipettor	0
	Washer	0
	Chambers	0
Cuvettes	BGW	1
	LC-le	1
	LC-ri	1

The requirement of system test results: details see in Chapter 9.1.

5. Daily Operating Process

This chapter mainly introduces the basic operating process of this system. After learning the content of this chapter, users are able to use this system to complete basic daily operation.

5.1 Daily Operating process

Table 5.1-1 Operating process for daily testing

Operating steps	Operation
1. Check before starting(1) Check the power supply(2) Check the analyzer	Confirm the power supply is normal; Make necessary check for the appearance and status of the analyzer.
2. Power-on	Connect the power supply , and turn on the main switch and submain switch of the analyzer. NOTE: auxiliary switch at analyzer MAGLUMI 800 is on left panel, MAGLUMI 2000 Plus has the same but on the front.
3. Log in the user software	Input the user name and password in [Login] dialog .
4. Check consumables	Confirm the system liquid, starter 1, starter 2 and cuvettes are enough to complete the test.
5. Confirm device status (1) Confirm incubator temperature (2) Confirm system test	Confirm whether incubator temperature is stable. The range of incubator temperature shall be $(36.8 \pm 0.5)^{\circ}\text{C}$ after it gets stable, with a fluctuation range within $\pm 1^{\circ}\text{C}$. Perform system test operation and confirm the tubing is full of system liquid and starter. Confirm whether the measured values of LC test and BGW test are within the allowable range.
6. Confirm test condition (1) Confirm assay (2) Profile/group setting	Confirm assay parameters. Set profile/group according to actual situation.
7. Prepare reagent (1) Confirm the remaining tests of reagent (2) Mixing time of magnetic microbeads	Confirm the remaining reagent tests of each assay is enough to complete the test. Confirm the mixing time of the magnetic microbeads is at least 30 minutes.
8. Assay calibration QC assay registration	Confirm whether the reagent has been calibrated, and perform calibration for the reagent without valid calibration. Perform QC operation for the reagent.

9. Sample registration	Perform conventional sample registration; Perform STAT sample registration; Perform dilution sample registration.
10. Test start	Click <Start> button to perform test.
11. Test of additional samples	Users can register additional samples during the test. (If the additional samples are STAT samples, STAT sample registration should be performed; if the additional samples are dilution samples, dilution sample registration should be performed); Click <Start> button to perform test.
12. Confirm test results	Search, recalculate, delete and print test results.
13. Operation at end of the test Quit operating software	Quit the software
14. Shutdown	Turn off the main switch and submain switch and disconnect the power supply of the device and the computer.
15. Operation after the test is finished	Take away the reagent from reagent area; Empty the waste tank; Remove the discarded cuvettes.

5.2 Test preparation

Complete the preparation work of conventional sample test via the following steps.

5.2.1 Check before Starting

Before starting, the following checks should be made to ensure the system works properly after starting.


NOTE

Make sure to wear gloves and work clothes to prevent infection. If necessary, wear protective goggles.

- 1) Check the power supply to confirm it supplies power normally. Check the communication cables and power cables of the analyzer, the computer and the printer to confirm they are connected properly.
- 2) Check whether the pipetting needle is at the correct position, whether the needle tip is dry, clean and smooth.

5.2.2 Power-on and Login the Software

- 1) Connect the power supply , and turn on the main switch and submain switch of the analyzer.

For MAGLUMI 2000 Plus

For MAGLUMI 800

***NOTE:** submain switch of MAGLUMI 800 is located on the left panel, for MAGLUMI 2000 Plus it is located on the front panel.

- 2) Connect the power supplies of computer and printer.
- 3) After logging in Windows operating system, double click the shortcut icon of the user software on the desktop. After starting, [login] dialog will appear on the screen. Enter the user name and password, and click <OK> button to enter the user software interface.
- 4) After the analyzer and all components complete initialization, wait until the incubator temperature gets stable. Then the test can be started.

NOTE

Maximum permissible time for setting the operating mode, calculated from switching on or starting the analyzer – 30 minutes

NOTE

The user name of system administrator is "snibe". Its initial password is "snibe".


Figure 2-1. [Login] dialog

5.2.3 Check Consumables

- 1) Check whether pipes for system liquid, starter 1 and starter 2 are connected properly and the liquid volume is enough;
- 2) Check whether there are enough cuvettes to complete the test;
- 3) Check whether the waste liquid is drained;
- 4) Check whether the discarded cuvettes are removed.

5.2.4 Analyzer status check

Confirm whether the device is normal via the following steps:

1. Confirm incubator temperature

Observe the color of button in the status bar. Test can only be started when the color is stable green. If incubator temperature exceeds $(36.8 \pm 0.5)^\circ\text{C}$ in the test process, the button turns red, but the analyzer will continue testing.

NOTE

After the analyzer is started, it will take about 30 minutes for the temperature in the incubator to be stable at $(36.8 \pm 0.5)^\circ\text{C}$. Therefore, wait for 30 minutes after startup of analyzer to perform sample test.

2. Confirm system test

Put the light check liquid on the rack with its tag facing the code reader, and load it to Channel 11 or 12 of the sample area (any channel for models 800). The light check liquid is automatically identified and displayed in yellow color in the Sample Info of the software, i.e. \$1c\$, as shown below.

Light Check inspection should be made once every week, or after any analyzer component is replaced. Do not load light check liquid into the analyzer unless the system performs the appropriate test.


Figure 2-2. [Sample area] interface

Click <System Test> button in the menu bar to enter [System Test] dialog. Here you can perform parameter setting for system test.


Figure2-3. [System test] dialogue

Input the priming times of the pipettor, the washer and the chamber. Input test times of BGW, LC-le and LC-ri.

Click **<OK>** button, [Message] dialog appears to confirm parameter setting, click **<OK>** button again.

Click **<Cancel>** button to cancel system test and exit [System Test] dialog.


Figure 2-4. Message

- **Priming:** Ensure the tubing of pipetting needle, washer and chamber is full of liquid.
- **BGW:** The reference range of BGW results is 200-1200..
- **LC:** The reference range of LC results is 400000-650000. , $CV \leq 3\%$. Difference between $LC(ri)$ and $LC(le) \leq 10\%$. (2000 Plus).

For specific expected values, refer to the instructions of the Light Check.

Test can be conducted only when BGW results and LC results meet the requirements.

In case of measurement abnormality, the test results might be not reliable.

5.2.5 Confirm Test Condition

Before performing tests, confirm that the reagent parameters and settings of this assay are correct. Perform the settings of group/profile.

1. Confirm assay parameters

Before using a new reagent, please read the instruction for use of the reagent carefully, and set up or confirm its parameters item by item according to actual requirements.

The set-up steps are as follows::

- Click **<Settings>** button in the menu bar, then click **<Test>** button.
- Select the assay to be edited, click **<Edit>** button to enter [User Specific Assay Data] dialog for test setting.
- Finally, click **<Save>** button to complete test setting.

2. Profile/group setting

In order to facilitate sample registration, the assays can be set with profile/group.

The set-up steps are as follows:

- Click **<Definition>** button in the menu bar, then click **<Profile>** button to enter [Profile Definition] dialog, where you can set the profile according to specific needs.
- Click **<Definition>** button in the menu bar, then click **<Group>** button to enter [Assay Group Definition] dialog, where you can set the group according to specific needs.
- .

5.2.6 Prepare Reagent

Scan reagent in the reagent area. Reagent calibration, QC test and sample test can be started if the reagent expiration date and the remaining reagent tests are satisfactory to perform the analysis and the mixing time of the magnetic reagent microbeads is at least 30 min.

1. Notes on use of reagent

ATTENTION!

Please remove the seals on the reagent kit before using a new reagent/

The preparation, use and storage of the reagent must be in strict accordance with its user manual. Prevent bubbles from forming in the reagent; or else the pipetting accuracy will be affected, hence affecting the test results.

The reagents of different kits cannot be mixed, or else the reliability of the test results might be affected.

2. Confirm remaining tests of reagent

After scanning RFID tag of the reagent in the reagent area and inserting the kit, make sure the remaining number of tests is enough to complete the number of tests needed. If it is not enough to complete the testing, please insert another reagent in the reagent area.

3. Shaking time of magnetic microbeads

When inserting the kit, the assay name, the remaining tests, the calibration information and the shaking time of magnetic microbeads are displayed in the window corresponding to associated buttons.

Figure 2-5. Reagent button

NOTE

The shaking time of the magnetic microbeads must be 30 minutes; or else the reliability of the test results cannot be ensured!

5.3 Test Analysis

After completing the preparation, the sample test can be performed.

5.3.1 Assay Calibration

Click <Reagents> button in the menu bar or in the [Home] interface, or open the reagent area door to enter [Reagents] interface.


Figure 3-1. [Reagents] interface

Confirm whether the reagent used in the test is calibrated and valid. If the reagent is without valid calibration, click **<Calibration>** or **<Batch Calibration>** button to execute calibration operation.

Valid calibration is required for the assay to calculate sample concentration.

After calibration, click **<View>** button to open **[Calibration Dialog]** dialog, where you can view working curve information and accept/reject new working curve.


Figure 3-2. [Calibration Dialog]

Select **Validate** and click **<OK>** button to accept the new calibration data.

Select **Reject** and click **<OK>** button to give up the new calibration data. Click **<Calibration>** or **<Batch Calibration>** button in **[Reagents]** interface to re-execute calibration operation.

Control Registration

1. Click **<Definition>** button in the menu bar, then click **<Control>** button in **[Definition]** menu to set up controls.


Figure 3-3. [Control] interface

2. Click <Add> button to open [Control Data Input] dialog.


Figure 3-4. [Control Data Input] dialog

3. Input the name and lot No. of the Control product, select the assay in need of control test, and click <Add> button to open [Control Detail Description] dialog.


Figure 3-5. [Control Detail Description] dialog

4. Input the expected concentration range, target value, target SD and replication times of this assay.
5. Click <Save> button, [Message] dialog appears to confirm parameter setting, click <OK> button to complete the control parameter setting.
6. Enter [Sample] interface. Properly load the rack with the control product, select the position of the control, and click <Control> button in the loading information area to open [Controls Selection] dialog, select the control information which has been set up, click <OK> button, and select "Start in Next Run" and click <Save> button to complete QC assay registration.


Figure 3-6. [Sample area] interface

**NOTE**

Sample test can be conducted only when the control test results of the reagent meet the requirements; or else the reliability of the sample test results cannot be ensured!

5.3.2 Sample Registration

**NOTE**

- 1) The samples with hemolysis, lipemia and icterus will affect the test results;
- 2) Ensure the sample is exclusive of clot; or else the pipetting needle will be blocked, which seriously affects the test results;
- 3) Some substances in the sample, such as drug, anticoagulant and preservative, will interfere with the test results;
- 4) Do not leave the sample open for a long time; or else the sample will volatilize, which affects the test results;
- 5) Improper parameter setting will affect the test results;
- 6) SNIBE Co., Ltd. recommends that the test results shall not be modified in any manner, and shall not be liable for any consequences arising there from.

Click <Sample area> button in the menu bar or in the [Home] interface, or open the sample area door to enter [Sample area] interface. Users can perform sample registration and view rack status.

Put the sample on the rack with its barcode facing the code reader, and load it to the samples area. The sample is automatically identified in the [Sample area] interface.


Figure 3-7. [Sample area] interface

Select the needed assay in the **Assay Selection** or select the assay in the **Profile Selection**. For STAT sample registration, select <STAT> button. The currently-registered STAT samples will be prioritized.

For dilute samples, select <Dilute> button, and select dilution ratio for the assay need to be diluted.

If the samples of the same rack are given the same assay, the whole rack can be selected so that the assay information can be added at one time for faster sample editing.

LIS application function can be used to obtain the assay info of the sample from the LIS server.

Click <Save> button to complete sample registration.

5.3.3 Start Test

After completing sample registration, click <Start> button in the menu bar to send the test command, and the analyzer will start test.

5.3.4 Additional Samples and Assays

1. Additional assay for registered samples

If additional assays are required for registered samples in the test process, assay registration should be first performed. In the **Rack Station** of the [Sample area] interface, select the samples that need additional assay, select the assay in the Assay Selection area, and click <Save> button to complete registration.

After the additional assay is registered, if there is no reagent for this assay in the reagent area, the reagent should be loaded; if there is reagent in the reagent area, click <Start> button in the menu bar to complete the continuous loading for the additional assay

2. Test of additional samples

If additional samples need to be added in the test process, sample registration should be firstly performed (refer to 5.3.3 Sample registration). Use a blank rack for the new samples.

Load the samples properly, select the needed assay in the **Assay Selection** area, and click **<Save>** button to complete sample registration.

For STAT sample registration, select **<STAT>** button. The currently-registered STAT samples will be prioritized.

To dilute samples, select **<Dilute>** button, and select dilution ratio for the assay need to be diluted.

If there is no reagent for this assay in the reagent compartment, the reagent should be loaded in the reagent compartment; if there is reagent in the reagent compartment, click **<Start>** button in the menu bar to complete continuous loading.

5.4 Test Results

5.4.1 Sample results

After the test, the test results can be searched in [Report] function. Processing of test results includes confirmation, printing, etc.

Click **<Report>** in the menu bar, then click **<LogI>** button on the left. Users can view, delete and modify the test results and print the journal.


Figure 4-1. [Journal] interface

Click **<Sort>** button and input corresponding date to search for historical journal information. User can also enquire journal information of certain ID, and choose to display test results according to the selected sorting criterion.

Click **<Valid>** button to confirm the selected test results according to the selected conditions. The confirmed journal is displayed in [Valid] interface.

].

5.4.2 Results and QC chart

After validation of the QC results, click the [QC] button on the [View Results] menu to go to the [QC] interface. Check the package QC chart or monthly QC process chart for the corresponding result as shown in the following figure.


Figure 4-2. [QC] interface

Select the "Package QC" or "Monthly QC" tab to go to the corresponding interface. Select the "Batch QC" or "Monthly QC" tab to enter the corresponding interface. User can select the time (time period for batch QC and specific month for monthly QC process), QC items, QC batch number to select relevant QC result data, and view the QC graph and QC analysis result generated by the program.

5.5 End of Analysis

When all the tests have been completed, quit the operating software and Windows operating system. Turn off the power supply of all parts.

5.5.1 Shutdown

When all the tests have been completed, quit the user software and Windows operating system. Turn off the power of all parts in the following order:

- 1) Turn off the power of the printer;
- 2) Turn off the power of the computer;
- 3) Turn off the power of the submain switch of the device.

NOTE

After the submain switch of the device is turned off, the reagent refrigerating system still works. Shutdown reagent cooling system should turn off the main power switch..

5.5.2 Operation after System Shutdown

To make preparation for the next test, perform the following:

- 1) Remove the samples in the sample area;
- 2) Remove the reagents in the reagent area;
- 3) Empty the waste tank;
- 4) Empty the discarded cuvettes in the waste bag;
- 5) Inspect whether the analyzer surface has stains. If so, wipe off the stains with a clean and soft cloth.

5.6[System] Menu Description

Click <System> button in the menu bar to enter [System] menu, where you can set up a series of functions, as described below:

Figure 66-1. [System] menu

Button	Function
<Info>	Display PC software, PLC software version info and analyzer serial number.
<Settings>	It consists of 4 tabs: "Mode", "Online", "User" and "Language". "Mode" tab: selection of the mode of operation and the mode of editing the sample information. "Online" tab: setting of parameters for connecting the software to the LIS server of the medical institution. "User" tab: changing the user name, password and access rights. "Language" tab: selection of the software interface language.
<Daily maintenance>	Set a daily maintenance reminder and display the priority of daily maintenance. Установка напоминания о ежедневном обслуживании и отображение порядка выполнения ежедневного обслуживания.
<Weekly maintenance>	Set a weekly maintenance reminder and display the priority of weekly maintenance.
<Monthly maintenance>	Set a monthly maintenance reminder and display the priority of monthly maintenance.
<Maintenance>	General cleaning and cleaning of the pipe system.
<Statistics>	Compiling statistics on reagent consumption.

5.6.1 <Mode>

Several sample editing modes and processing modes are available to the user in this software. The settings can be changed according to the actual conditions.

Click the Mode tab to enter the [Mode] interface as shown in the figure below.


Figure 5.6-2. [Mode] interface

[Mode] interface displays the following information:

1. Running Mode

- Random Access Mode: processing one rack after another, prioritized as follows:
 - A. STAT assay;
 - B. Assays with automatic reference and automatic dilution;
 - C. From left to right racks;
 - D. According to incubation time, from the longest to the shortest;
 - E. According to assay name abbreviations, from A to Z;
 - F. According to positions of sample pipes on racks;Advantage: this mode allows removal of a rack after the whole rack is process and re-adding of new samples.

Disadvantage: processing time is not optimal

- Batch mode: Processing all samples with optimal time, prioritized as follows:
 - A. STAT assay;
 - B. Assays with automatic reference and automatic dilution;
 - C. According to incubation time, from the longest to the shortest;
 - D. According to assay name abbreviations, from A to Z
 - E. From left to right racks, from back pipes to front pipes within a rack.Advantage: This mode makes full use of cuvettes with the highest efficiency.
- Disadvantage: Unable to process a whole rack so as to add new samples.

2. Edit Sample Mode:

- Normal Mode:
Barcode reader directly reads the barcode from sample pipes to generate sample ID.
- Quick Mode:
Software directly generates sample ID without reading barcode from sample pipes.
- Online Mode:
Software reads assay information from Lis server via sample ID, and sends test results to hospital Lis server.
- Emergency Mode:
When there is an error with the sample area door sensor or barcode reader (after the rack is loaded, the corresponding position of “rack station” in [Sample area] interface displays “ERROR”) and thus “sample info” cannot be edited in the sample loading interface, emergency mode can be enabled. In this mode, users can still edit sample info.

5.7 [Settings] Menu

The settings for performing basic tests are set by the manufacturer. These settings can be changed using the function keys in the [Settings] menu. As a rule, the set parameters apply to all tests.

Select <Settings> from the main menu to go to the [Settings] menu. The function buttons are described below.


Figure 5.7-1. [Settings] menu

Button	Function
<Test>	Detailed configuration of assay parameters.
<Control>	Controls setting.
<Group>	Control group setting.
<Class>	Defining multiple assays as a class. When performing one type of assay on multiple samples, they can be combined into a class for aggregate processing.
<Delutant>	Setting the assay in which it is required to perform dilution and specify the dilution factor.
<Results setting>	Settings for displaying, editing, and printing results.

5.8 [Process] menu

Select <Process> menu in the menu bar so to move to [Process] menu, where a number of fuctions described below shall be set.


Figure 5.8-1. [Process] menu

Button	Function
<Initialization>	Initialize the analyzer.
<Init W. Clear>	Initialization with Cuvette(s) Clear/
<Restart >	Restart uncompleted assays/
<Return Asy>	Return uncompleted test time.
<Warning Opt.>	Set up the warning message in emergency circumstances

5.9 [System test] menu

[System Test] menu is used to wash pipe system of the analyzer. BGW and LC are used to detect pipetting system, washer, chamber and starter. All system test must be performed before the analyzer starts an assay. Click [System Test] button in the menu bar to enter [System test] dialog


Figure 5.9-1. [System test] dialog

[System Test] dialog contains:

- **Pipettor:** Set up wash cycles for the pipe, hose and pipetting needle of the pipetting system.

- **Washer:** Set up wash cycles for the wash pump, hose and wash needle of the washer. Users must perform wash before the analyzer runs.
- **Chamber:** Set up wash cycles for the pipe system of the chamber.

2. Cuvette

[BGW]	Enter BGW cycles (default value: 1).
[LC-le]	Enter cycles to test the left pipettor light signal. Use the left pipettor to dispense liquid to test LC (default value: 1).
[LC-ri]	Enter cycles to test the right pipettor light signal. Use the right pipettor to dispense liquid to test LC (default value: 1).

BGW and LC must be performed after any starter is replaced.

The requirement of system test results is:

- 1) BGW: RLU is 100-1200, $CV \leq 10\%$;
- 2) LC: RLU is 450000-650000, $CV \leq 3\%$;
- 3) For two-pipettor modes: mean difference between LC(ri) and LC(le) $\leq 10\%$.

Please refer to the LC Reagent Manual of the specific expected value.

To perform BGW and LC, the default cycle (1) can be changed according to needs. Meanwhile, the cuvettes of corresponding number should be put in. If BGW or LC is not necessary, the default cycle (1) can be changed to (0).

Click **<OK>** button, [Message] dialog appears to confirm parameter setting, click **<OK>** button again

Click **<Cancel>** button to cancel system test and exit [System Test] dialog.


Figure 5.9-2. [Message] Dialog for System Test

NOTE

1. Complete system check before the test.
2. When selecting LC-le / LC-ri, the rack carrying the LC liquid must be placed on Track 11 or 12 on the right in the sample area! (without limitation for MAGLUMI 800)

When the analyzer is running an assay, if you want to enter **<System Test>** from the menu bar, a message confirmation dialog appears to prompt that system test can only be performed after end of the assay.


Figure 5.9-3. [Message] Dialog to cancel operations

5.9.1 Load Light Check

Light Check produced by SNIBE Co., Ltd. is used for testing the pipetting needle performance of analyzer. When performing Light Check, put prepared Light Check liquid tank on the rack with its tag surface facing the code reader and load it to channel 11 or 12 in the sample area. The liquid will be automatically identified and displayed in **Sample Info** of the software in yellow, i.e. \$!c\$.

Manually run LC sample

If the code reader fails or the tag of Light Check liquid bottle is damaged, manual editing can be made in the software. Put LC liquid bottle on the tube rack and load it to channel 11 or 12 of the sample area (without limitation for MAGLUMI 800). Use the mouse to select the position of Light Check liquid in the sample area; Light Check must be defined by pressing **<LC>**, with this position displayed in yellow, i.e. **\$LC\$**. Click **<Save>** to return to the **[Home]** interface, and click **<Start>** to perform LC check.

The requirement of LC check results is:

- 1) LC: RLU is 450000-650000, CV $\leq 3\%$;
- 2) For two-pipettor modes: mean difference between LC(ri) and LC(le) $\leq 10\%$.

Please refer to the LC Reagent Manual of the specific expected value.

5.9.2Add Light Check Liquid

Light Check liquid is used for LC check.

LC should be performed under the following three circumstances:

- 1) Before start of the first round of test on a daily basis;
- 2) After change of the Lot of starters;
- 3) After maintenance of the analyzer.

Prepare Light Check solution:

1. Take a bottle of Light Check out of the packaging box.
2. Carefully unplug the bottle.
3. Prior to use, the Light Check solution must be placed still for at least 5min. Light Check test should be performed at room temperature.
4. Insert the Light Check bottle (with rubber plug removed) to the sample rack. The bottle barcode must be readable. Then insert the sample rack into channel 11 or 12 in the sample area of the analyzer.

WARNING!

Prevent air bubbles from forming, which may affect the Light Check result in

5.10 [Reagents] interface

Users can open the [Reagents] menu in the following 3 ways:

- Click **<Reagents>** on the menu bar;
- Open the reagent area door;
- Click the **Reagents** icon on the **[Home]** interface.

After entering the interface, users can:

1. Load or unload reagents;
2. View/ accept/reject/validate calibrator results and start calibrate


Figure 3.10-1. [Home] interface

On the [Home] interface, Reagents icon is mainly shown in five colors:

- **Red:** Not recognized by the system;
- **Yellow:** Recognized by the system, without valid calibration data;
- **Green:** Recognized by the system, with valid calibration data;
- **Purple:** Recognized by the system, with expired calibration data;
- **Black:** Recognized by the system, but with expired reagents.

Click **<Reagent>** button in the menu bar to enter **[Reagents]** interface.


Figure 5.10-2. [Reagents] interface

Reagent Area	Shows the status and name of the reagent.
Reagent Data	Shows the relevant data of the reagent cartridge.
Calibration	View the calculate curve of the reagent. Start calibrate.
Reagent Operations	Exit [Reagents] interface. Perform batch calibration.

5.10.1 Reagent Area

Reagent Area shows the status and name of the reagent in the reagent area, the magnetic microbeads mixing time, etc.

The digits below the assay name indicate the number of tests that can be performed with the reagent. The countdown timer indicates the magnetic microbeads mixing time.

1. Definitions of numbers and colors

Red and green are used to show the state of reagents and give prompts as to when the reagent can be taken out of the reagent area.

Table 5.10-3. Reagent status

	Number 1 means the reagent is located in No. 1 track of reagent area. Green means the reagent is not used and can be taken out of the reagent area.
	Number 2 means the reagent is located in No. 2 track of reagent area. Red means the reagent is being used and cartridge cannot be taken out of the reagent area.

NOTE


Users must pay attention to the usage state of each reagent!!

2. Reagent loading conditions

There are four types of reagent loading conditions.


Table 5.1010-4. Reagent loading Conditions

	No reagent inserted into this channel.
---	---


	Система распознала реагент, он находится в канале № 2; цвет фона светло-серый, что значит, реагент не выбран.
	The reagent has been recognized by the system; the reagent is located in No. 2 track of reagent area; the background color is dark gray, indicating the reagent is selected.
	Reagent is available in the channel, but has not been recognized by the system. The reagent needs to be reloaded.

3. Status of calibration

Two small squares with different color display the status of calibration.

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Status of Calibration displayed on the screen for MAGLUMI 2000 Plus:

 ООО ЦИРМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU

Table 5.10-5. Calibration status

Symbol		Definition
 (Red)	 (Grey)	Without valid calibration.
 (Red)	 (Red)	Without valid calibration, calibration is in progress.
 (Red)	 (Green)	Without valid calibration, new calibration has been executed but has not been validated.
 (Green)	 (Grey)	With valid calibration (calibration has taken effect).
 (Green)	 (Red)	With valid calibration, the second calibration is in progress.
 (Green)	 (Green)	With valid calibration, and the second calibration has been executed but has not taken effect.
 (Purple)	 (Grey)	The calibration has expired.
 (Purple)	 (Red)	The calibration has expired; the second calibration is in progress.
 (Purple)	 (Green)	The calibration has expired; the second calibration has been executed but has not been validated.
 (Black)	 (Black)	The reagent has expired.

4. Status of time

The countdown timer indicates the magnetic microbeads mixing time, which is 30min by default. If you click <Start> before the countdown timer reaches 00:00, the analyzer will give a prompt message indicating that magnetic microbeads are not homogeneously mixed. In such case click the <OK> button to start test.

The magnetic microbeads mixing time will be reset to 30:00 in the following two ways:

- 1) after the analyzer is restarted.
- 2) the reagent is reloaded after it has been taken out for more than 2min.

5. Order of using reagents

When there are more than one reagent of the same type loaded in the reagent area, such reagents should be used in the following sequence:

- 1) Priority is given to reagents allowing fewer remaining tests.
- 2) Use reagents from left to right.

5.11 [Reagent data] menu

[Reagent data] interface can be entered in three ways:

- open the door in the sample area;
- click the Sample area icon on [Home];
- click <Sample area> in menu bar.

Possible user actions after entering the “patients” interface:

- 1) Loading the following sample types for analysis: sample, control, and LC liquid;
- 2) Selecting the assay for the sample.
- 3) Setting STAT test or dilution test.
- 4) Editing the sample details.


Figure 5.11-1. [Home] interface

The “patients” histogram can be displayed in three colors:

- Red: the installed rack is not recognized by the system.
- Green: the installed rack is recognized by the system.
- Grey: no racks installed.

[Sample area] interface consists of the following elements:


Figure 5.11-2. [Sample area] interface

Rack station	Display of rack loading status.
Loading	Specifies the type of test, such as “STAT”, “Control”, “Dilution Selection”, or “LC check liquid”.
Sample data	Displays sample number.
Select test	Allows users to select assays for different samples.
Select class	Allows users to quickly find the sample class they are looking for.
<Save>	Exit [Sample area].

In the “Select test” area options BGW, LC-le and LC-ri are used for system.

5.11.1 Rack station

After a rack is added, relevant sample data will be read by the barcode reader. Each button represents a sample channel. Select a rack number, and press the button; the sample ID will be shown in the **Sample Info** area, and the sample channel is shown in blue.


Figure 5.11-3. Rack station

The meanings of different colors are listed below:.

	This icon indicates the sample channel is empty.
	This icon indicates the channel has valid samples and has been selected.
	This icon indicates the channel has valid samples but has not been selected.
	This icon indicates the rack inserted in the channel is not recognized and must be reinserted.

5.12 [Report] menu

Click **<Report>** on the menu bar to open **[Report]** menu. The buttons are on the left:


Figure 5.12-1. [Report] menu

Button	Function
< Journal>	display the journals on the then-current day, and search and display historical journals.
<Valid results>	display the validated journals on the then-current day, and search and display historical valid results.
<Calibrator>	display the validated calibrator results on the then-current day, and search and display historical calibrator results.
<QC results>	display the validated control test results on the then-current day.
<System test>	display all system test results.
<QC>	display content related to quality control chart.
<Report>	display the assay report and relevant settings.

6 Status Display and Handling

6.1 Status display

Figure 1-1. Status bar

Software interface at the bottom of the status bar has three buttons to display the status of the analyzer. They are **<Reservoir>** () , **<Temp/Volt.>** () and **<Waste>** () .

6.1.1 Reservoir

The **<Reservoir>** button is used to estimate the number of cuvettes, wash liquid level, and starters 1 and 2.

The colors of the **<Reservoir>** button indicate the following:

green

All consumables are sufficient.

Yellow

One or more consumables are not sufficient.

Red

One or more consumables are almost used up.

WARNING!

In case of an alarm, please check which consumable has run out and replenish it.

Click **<Reservoir>** button to open [**Reservoir Status**] dialog, where shows the status of cuvettes , starter reagents and system liquid.


Figure 1-2. [Reservoir status] dialog

6.1.1.1 Cuvettes

Green	MAGLUMI 800 Number of cuvettes bars: >5. MAGLUMI 2000 Plus Number of cuvettes bars: >17.
Yellow	MAGLUMI 800 Number of cuvettes bars: from 1 to 5 inclusive. At this moment, a warning prompt "Running out of cuvettes soon!" will appear, and an alarm sound will be given. MAGLUMI 2000 Plus Number of cuvettes bars: from 4 to 17 inclusive. At this moment, a warning prompt "Running out of cuvettes soon!" will appear, and an alarm sound will be given.
Red	MAGLUMI 800 Number of cuvettes bars: 0. At this moment, a warning prompt " Run out of cuvettes! " will appear, and an alarm sound will be given. The analyzer will report an error and shut down. MAGLUMI 2000 Plus Количество пластин кювет: < 4. At this moment, a warning prompt " Run out of cuvettes! " will appear, and an alarm sound will be given. The analyzer will report an error and shut down.

6.1.1.2 Wash liquid

Different volume of wash liquid show by following icon.

Normal	Volume of wash liquid: ≥ 20 %.
Almost empty	Volume of wash liquid: <20 %.
EMPTY	Wash liquid has run out; the analyzer reports an error and shuts down.

WARNING!

1. Turn off the analyzer before adding wash liquid. The user can connect the wash liquid bottle to another bottle using the continuous feed tube to replenish the liquid level while the analyzer is running.
2. For information on storing and using wash liquid, see the

instructions in the wash liquid package.

6.1.1.3 Starters

Different volume of starters show by following icon.

Normal	Volume of starter 1 and 2: $\geq 20\%$.
Almost empty	Volume of starter 1 and 2: $< 20\%$.
Empty	Starter 1 and 2 has run out; the analyzer reports an error and shuts down.

WARNING!

The starter 1 and 2 can not be exchanged when the analyzer is testing.

6.1.2 Temperature and Voltage

The status of <Temp/Volt> color button means:

All temperatures and voltages in the analyzer are normal.

One or several temperatures or voltages in the analyzer are beyond the normal/

Click <Temp/Volt> button to open [System Parameter] dialog, which shows the temperature and voltage parameters of the analyzer.


Figure 1-3. [System parameter] dialog

6.1.3 Waste

The status of <Waste> color button:

The status of all kinds of waste are normal.

One or more kinds of waste is nearly full.

One or more kinds of waste are full.

Click **<Waste>** button to open **[Waste Status]** dialog, where shows the volume of liquid in the waste liquid tank and the number of discarded cuvettes in the waste bag



Figure 1-4. **[Waste tank status]** dialog

6.1.3.1 Liquid waste

Different volume of waste liquid show by following icon.

Volume of liquid in the waste liquid tank ≤ 80 %.

Volume of liquid in the waste liquid tank reaches 80–90 %.

Volume of liquid in the waste liquid tank reaches 90–100 %.

WARNING!

When yellow or red appears, the system will display a warning message; in such case, please empty the waste liquid tank.

6.1.3.2 Reaction module

Green	MAGLUMI 800
	Number of cuvettes bars discarded: <60 .
Yellow	MAGLUMI 2000 Plus
	Number of cuvettes bars discarded: <100 .
Yellow	MAGLUMI 800
	Number of cuvettes bars discarded: from 60 (inclusive) to 80 (not inclusive).
Yellow	MAGLUMI 2000 Plus
	Number of cuvettes bars discarded from 100 (inclusive) to 120 (not inclusive).
Red	MAGLUMI 800
	Number of cuvettes bars discarded: ≥ 80 .
	MAGLUMI 2000 Plus
	Number of cuvettes bars discarded: ≥ 120 .

**WARNING!**

When yellow or red appears, the system will display a warning message; in such case, please empty the waste liquid tank.

6.2 Reservoir Handling

In cases where consumables are running low or waste containers are filling, the analyzer will display a warning message and shut down, and the [Operation Suspended] dialog will appear in the software interface. Once the consumables level has been recovered and waste has been removed, click the <Continue> button to resume assay (initialization of stacker is required).

6.2.1 Placing Cuvette

Operation Steps:

1. Unpack the cuvettes as instructed, and take out a layer of (8) cuvettes bars.
 2. Place the cuvettes on the conveyor belt along the proper direction.
 3. The conveyor belt will automatically transfer the cuvettes into the stacker.
 4. When the conveyor belt stops, users can continue adding cuvettes.
 5. Repeat the above steps until the stacker becomes full.
- Cuvettes can be added at any time when the analyzer is running.

**NOTE**

- 1) Pay attention to the remaining number of cuvettes during assay!
- 2) The stacker should be emptied once a week to ensure cleanness of cuvettes.

6.2.2 Adding Wash Liquid

Wash liquid is connected at the right side of the analyzer. "Wash liquid" is used for cleaning pipette, system pipe and magnetic microbeads. The volume of wash liquid is detected by the level detector, and shown with the "Reservoir" on the main menu. Click <Reservoir> button to open [ReservoirStatus] window. The displayed information includes the stock of cuvettes in the stacker, starters and wash liquid.

**WARNING!**

Make sure there is sufficient wash liquid before and during the assay!

The way of add System Liquid in standby:

1. Take the sealing cover with level detector and hose out of the corresponding bottle.
2. Add wash liquid, and close the sealing cover.
3. Fill the hoses: Select <System Test> on the menu bar, and select pipettor and washer in the [System Test] dialog. Enter (3) and (6) in corresponding cycles, and click <OK> to start priming.

The way of add Wash Liquid in testing:

Use the Continuous Loading Pipe and another System Tank from the Wash Liquid kit to add system liquid during the test.

6.2.3 Change Starter Reagents

The container of starter is on the right side of the analyzer. When you open the enclosure door, you will see the S1 and S2 marks on the reagent pipes. The digit S1 and S2 represent starter 1 and starter 2, respectively.

The reagent bottle is sealed with an openable screw cap. The level of starter can be detected with the level detector fixed to the cap.

WARNING!

- 1) Make sure there is sufficient starter before operation of the analyzer!
- 2) The starters can not be changed when the analyzer is testing.
- 3) Please make sure hoses for starters 1 and 2 are connected correctly.

WARNING!

- 1) Never pour out starters, liquid splash is not allowed in this area!
 - 2) Never mix starter 1 and starter 2 manually.
 - 3) Chemical burn risk! Please read the reagent instructions in the package of starter.
-

Operation steps:

Change starters of the same lot number:

1. Open the starter 1 bottle and take out the connecting hose.
2. Place a new bottle of starter 1 at the corresponding position.
3. Close the screw cap.
4. Replace starter 2 according to the above steps.
5. The connection pipes of the starter are marked with different numbers. Make sure that the pipes are connected correctly.
6. Fill the hoses: Select **<System Test>** on the menu bar and select chamber in the [System Test] dialog. Enter (3) in corresponding cycles, and click **<OK>** to start priming.

Change starters of different lot number, please do the following additional operations:

1. Background priming: Select **<System Test>** on the menu bar. Change (1) preset in the BGW, Lc-le and Lc-ri area to (0), and press **<OK>** to start the BGW test.

6.2.4 Use of cuvette

Waste cuvette holder is positioned on the right side of the analyzer. The software will automatically count the number of waste cuvettes and prompt users.

When the waste bag is full, users can take it out of the holder. Replace with a new waste bag and fit back in the holder. Click **<Waste>** button in status bar to open [Waste Tank Status] dialog, click **<Reset>** button to reset the number of waste cuvettes.

WARNING!

- 1) After replacing the waste cuvettes, please reset the number of waste cuvettes or the warning message will always remain. Or the warning message will always remain displayed.
-

-
- 2) Make sure the waste bag is placed properly; otherwise the analyzer may be disrupted as the cuvettes are stuck at the edge of the waste bag when they are pushed out of the chamber.
 - 3) The waste cuvettes shall be placed in the waste bag since residual substances in the cuvettes may cause contamination. Please dispose waste cuvettes in accordance with local laws and regulations.
-

6.2.5 Liquid waste

The liquid waste port is located inside the analyzer enclosure. If you need to remove liquid waste, open the front enclosure door, remove the tray, unscrew the red liquid waste tank lid, remove it, and remove the liquid waste.

Waste liquid comes from cuvettes, pipetting system and washer during assay, including magnetic microbeads, starters, samples, reagents, and wash liquids.

WARNING!

- 1) Please wear gloves for operation!
 - 2) Please dispose waste liquid in accordance with local laws and regulations.
 - 3) Please clean the waste liquid tank on a regular basis according to the Maintenance Manual.
-

7 Reagent Loading

7.1 Reagent cartridge

All reagent cartridges basically have the following structure:


Figure. 1-1. Reagent cartridge

1. Seal is for sealing reagent. The pipetting needle enters the kit from here to absorb reagent;
2. Handle is convenient for users to hold the kit by hand;
3. The RFID label which is behind the ID label records the reagent data;
4. The ID label indicates the reagent components, name, expiry date and other info of the kit;
5. The magnetic microbeads mixing gear is used in conjunction with the rocker in the reagent area to mix magnetic microbeads homogeneously;
6. ID Tab is used to block optocoupler when the reagent is inserted to the reagent area so that it can be identified;
7. Clamp is used to fix the reagent inserted to the reagent area.

7.2 Reagent loading

7.2.1 Reagent preparation

It is not allowed to store reagent upside down. Besides, please do not shake the reagent.

1. Take the reagent out of the cooling chamber;
2. Take the reagent kit out of the package;
3. Check the sealing film and other parts of the reagent kit to see if there is any leakage. In case of leakage, please contact your local agent immediately;
4. Check each component liquid of the reagent kit to see if air bubbles exist. If yes, please use a pipettor to remove the air bubbles; the reagent integral can be used after confirming that all air bubbles are eliminated;
5. Gently turn the magnetic microbeads mixing gear to make sure it can rotate freely;
6. Carefully tear off the kit sealing film;
7. Clean up liquid on the reagent kit surface to avoid nosocomial infection.

7.2.2 Reagent loading

Before inserting the reagent to the reagent area, please make sure the software and the computer (analyzer) are working normally.

1. Open the reagent area door; the software will automatically display the [Reagents] interface;

2. Hold the reagent handle to get the RFID label close to the sensing area (for about 2s); the buzzer will beep; one beep sound indicates successful reading;
3. Keeping the reagent straight insert it along the blank reagent track;
4. Relevant reagent data will be read out via computer software; the [Reagents] interface will show the reagent name, and the corresponding reagent position button will change to dark gray;
5. In case of failure to read data, please repeat the above operations;
6. After finishing reading reagent data, please close the reagent area door and exit the [Reagents] interface.


Figure 2-1. Reading and loading the reagent

7.2.3 Reagents removal

When taking the cartridge out of the analyzer, make sure both the software and the analyzer are running normally, so as to ensure the reagent is removed in normal condition. The reagent should not be removed if there are still some assays in progress or pending for results on the analyzer (the corresponding channel in the [Reagents] interface will be highlighted in red). Otherwise the result will be an error, and may cause damage to the analyzer.

1. Open the reagent area door;
2. Hold the reagent handle to take the kit out of the reagent area;
3. Do not turn the cartridge with the reagents upside down.

If the reagent is empty

1. Please properly dispose the reagent.

If the reagent is not empty

1. Put the cartridge and store it in a cooling chamber, ensure that the up and down orientations of the reagent are correct.

7.3 Proper reagents storage

1. Keep the reagent with correct up and down orientations ensured;
2. If the reagent is opened, please use additional sealing film to cover each hole of the reagent in order to avoid liquid evaporation.

8 Sample Loading

8.1 Sample rack

The sample rack specified for the analyzer has 12 sample positions (for MAGLUMI 800 — 8 positions). The sample tube can be barcoded for identification. The program automatically reads the barcode.

8.1.1 Sample rack assembly

The sample rack consists of the following components:

1. Operating handle: is used for loading and unloading the sample rack;
2. The sample position barcode is used to identify the position of the tube; there are also fixing clamps on the handle, numbered 1 to 12 (1 to 8 for MAGLUMI 800 analyzers).
3. Sample tube position: has a retaining clip inside to fix a sample tube ;
4. Rack fixing clamp: can be used to fix the sample rack and detect the status of sample tube.

The barcodes on the sample rack are used to determine the specific positions of sample tubes. When sample rack is inserted, the barcode reader will scan tube info on the rack one by one/ The position of each sample and other info contained in the barcode will all be displayed on the Sample Info of [Sample area] interface.

After the sample rack is correctly inserted to the sample reservoir, the sample info will be read.


Figure 1-1. Sample Position info

No. 3 & 4 sample tubes have no barcode, therefore cannot be automatically identified. If a sample is not automatically identified, users can manually edit the sample info.

8.1.2 Sample rack type

The analyzer only uses one type of sample rack. Sample tubes with the inner diameter of $\varnothing 8$ mm– $\varnothing 12$ mm and the height of 60 mm–100mm can be used.

8.2 Loading Sample

WARNING!

- 1) Please make sure the sample tube has been uncapped before pipetting; otherwise the sampling needle will get damaged.
- 2) In order to avoid infection, please wear gloves before using the sample.

8.2.1 Sample preparation

NOTE

In case of lack or absence of sample solution, the analyzer will still run the test; If the system displays a warning message such as “Left/Right Needle can not detect Liquid!”, then you must run the test again.

The following conditions must be satisfied before measuring sample:

1. Samples include: serum, plasma, urine o whole blood, etc.;
2. Prepare anticoagulants such as EDTA and heparin;
3. The blood drawing process should be standardized; serum should be separated from blood clots;
4. If the sample contains suspended solids or is turbid, or has blood lipids or RBC fragments, filtering or centrifugation must be performed prior to test;
5. Samples with hemolysis or lipemia, or infected by microorganism, should not be used for test;
6. Please ensure the sample contains no air bubbles prior to measure.

NOTE

To ensure safety, the sample must satisfy relevant requirements and operation conditions to avoid air bubbles and blood clots!

The sample satisfying all conditions can be put in a sample tube to be inserted to the sample rack. The operation steps are described as follows:

1. Test tubes must conform to specification requirements;
2. Carefully insert each sample tube to the sample rack;
3. If barcode scanning will be performed for the sample rack, each barcode should face rightward so that it can be read by the barcode reader.

NOTE

- 1) Do not revolve a sample tube with barcode in the sample rack to avoid damaging the barcode.
- 2) After placing the test tube rack in the sample area, do not change the position of the tubes randomly: this will confuse the test results.

Like samples, each barcode of Light Check fluid tank provided by Snibe should face the barcode reader so that it can be read. If a sample tube is not properly inserted to the sample rack, please pull it out and reinsert it

8.2.2 Loading sample rack

Please load the sample rack according to the steps below:

1. Open the door of sample area;
2. Select an empty slot to push in the sample rack;
3. Slowly insert the sample rack so that each barcode can be read effectively;
4. Make sure the sample rack is kept vertical when it is pushed in; push in the rack smoothly until the edge of sample rack is in close contact with the edge of sample area;
5. In case of barcode scanning failure, please reload the sample rack and scan the barcode.

NOTE

- 1) Before pushing the sample rack into the sample compartment, please make sure both the analyzer and the software have entered the working status. Otherwise, the sample cannot be read and subsequent operation will be invalid.
- 2) The barcode reader has laser radiation hazard. Please avoid looking directly at the laser beam!

Remove the Sample Rack

Before removing the sample rack, please make sure both the analyzer and the software are running normally. The sample rack to be removed should not be in active status (orange light is on); otherwise, a result error will occur and the analyzer may get damaged.

1. Open the door of sample area;
2. Grasp the rack handle with your thumb and index finger to smoothly pull the sample rack out of the sample area until the sample rack completely leaves the sample area.

If the tube is empty

3. Dispose the sample tube according to relevant regulations.

If the sample tube is not empty and may be used in the future

4. Seal and store the sample properly according to relevant laboratory regulations.

8.3 Maintenance of test tube racks

The maintenance of sample rack is necessary under the following special circumstances.

1. The rack is in direct contact with any biochemical substance.
2. A sample tube cannot be inserted to the sample rack stably.

In the first case:

1. Transfer all samples on the contaminated sample rack to another clean rack;
2. Soak the contaminated sample rack in 0.1% sodium chloride solution for 15~30min. In order to avoid corrosion of the sample rack, the duration should not exceed 30min;
3. Take out the sample rack; use a clean, dry towel to wipe it dry.

In the second case:

1. Take a ball-point pen; insert its tip below the retainer of metal plate inside the sample rack slot, as shown in Fig. 9.3-1; then, apply gentle force to make the retainer upward;

2. Insert a sample tube conforming to specification requirement and check if the problem is fixed;
3. If the sample tube is still loose, please repeat Step 1 and 2 until the tube can be inserted firmly.


Figure 3-1. Adjust the retainer inside the slot

9 System Maintenance

9.1 Analyzer Maintenance

9.1.1 Daily Maintenance

Materials needed:

Clean waste bag, clean cotton cloth.

Maintenance Process:

1. Use the clean cotton cloth to clean the analyzer surface;
2. Empty the waste bag to remove waste cuvettes and replace for a clean bag;
3. Check the volume of starters; if not enough for use of the next day, please replenish or change it;
4. Check the volume of wash liquid; if not enough for use of the next day, prepare a new wash liquid.
5. Perform the "Priming for All" function once more time.

9.1.2 Weekly Maintenance

Materials needed:

Cotton swabs, alcohol.

Maintenance Process:

Use cotton swabs dipped in alcohol to clean the outer surface of pipettor, washer's injecting needle and aspirating needle; then perform the "Priming for All" function twice.

9.1.3 Monthly Maintenance

Materials needed:

Sodium hypochlorite solution with effective chlorine content of 0.05%~0.1% (i.e. 500 mg/L~1000 mg/L), distilled water, system tubing cleaning solution, light check liquid.

Maintenance process:

1. Use 40 °C distilled water to replace starters, and perfuse the chamber 30 times; then replace the water with starters, and perfuse the chamber 10 times;
2. Use System Tubing Cleaning Solution to perform the "Wash Pipe" function once;
3. Empty the Wash Liquid Tank and the waste liquid Tank:
 - 3.1 Empty the wash liquid tank and the liquid waste tank;
 - 3.2 Add 2 liters of sodium hypochlorite solution with an effective chlorine content of 0.05%~0.1% (i. e. 500 mg/L~1000 mg/L) to the drum, cover with a lid, shake back and forth, and leave for 30 min;
 - 3.3 Pour the cleaned sodium hypochlorite solution and add clean water (2 liters each time) to the cleaning bucket at least 3 times.
4. Fill a bucket with cleaning solution to spray the pipettor and washer station, 30 times each;
5. Place Light Check Liquid container in the upper right corner of the sample compartment. Perform the BGW test and the LC once each. Review and record the results.
6. System test results requirements:
 - (1) BGW: RLU – 100–1200, CV ≤10 %;
 - (2) LC check: RLU – 450000–650000, CV ≤3 %;

(3) Dual metering modes: average difference between LC-ri and LC-le $\leq 10\%$ (2000 Plus).

(4) Please refer to the LC Reagent Manual of the specific expected value..

9.1.4 Lasting System Shutdown

If you do not plan to use the analyzer for 5 or more day, the following steps should be completed:

1. Replace the Starter 1 and 2 bottles with the two bottles of distilled water.
2. Replace the system liquid tank with a tank of distilled water.
3. Click <**System Test**> button in menu bar, select Pipettor, Washer and Chamber, the cycles is 30, 20 and 10 respectively.

When using the analyzer again, the following steps should be completed:

1. Remove the distilled water bottles and the distilled water tank. Click <**System Test**> button in menu bar, select Pipettor, Washer and Chamber, the cycles is 30, 20 and 10 respectively.
2. Place two bottles of Starter 1 and 2.
3. Place a tank of wash liquid.
4. Click <**System Test**> button in menu bar, select Pipettor, Washer and Chamber, the cycles is 30, 20 and 10 respectively. System test results shall meet the requirements

9.1.5 List of Components

NO.	Component description	Location	Replacement period	Notes (component replacement)	Relevant models
1	Hose pump tube for pumping liquid waste out of the washer station	Top side of the washer station	Requires replacement once a year and as needed if damaged	Remove the pump and disconnect the inlet and outlet ends of the tube to be replaced. Then open the semi-circular clamping device for the hose pump tubing so that it can be removed. First remove the used tubing, then cut a new tubing of the same length and connect it to the pump so that the direction of the new tubing matches the direction of the removed tubing (the direction of the retainer). Tighten the retainer and insert the inlet and outlet ends of the tubing one after the other into the mating holes. To replace the remaining tubes, follow the sequence of operations described above. Finally, connect the hose pump to a power outlet.	MAGLUMI 800
2	Pipettor probe	Pipettor	Requires periodic replacement of components and replacement in case of bending damage	1. First remove the housing of the lower section of the pipettor, then the small plate of the communication cable and then disconnect the tube ends. After that, loosen the small metal stopper by turning it in reverse direction. Next, loosen the plastic nut located under the pipettor probe so that the probe module can be removed. Finally, loosen the small nut and remove the set	

System maintenance

				collar and spring to remove the pipettor. 2. Insert a new probe into the set collar and spring, tighten the small nut, insert the probe module into the pipettor and tighten the plastic nut (be careful when positioning the set collar). Rotate the small metal stop to lock the position of the probe, then secure it. Finally, insert the ends of the tubing, secure the communication cable and install the housing of the lower section of the pipettor.	
3	Hose pump tube for pumping liquid waste out of the washer station	Top side of the washer station	Requires replacement once a year and as needed if damaged	Remove the pump and disconnect the inlet and outlet ends of the tube to be replaced. Then open the semi-circular clamping device for the hose pump tubing so that it can be removed. First remove the used tubing, then cut a new tubing of the same length and connect it to the pump so that the direction of the new tubing matches the direction of the removed tubing (the direction of the retainer). Tighten the retainer and insert the inlet and outlet ends of the tubing one after the other into the mating holes. To replace the remaining tubes, follow the sequence of operations described above. Finally, connect the hose pump to a power outlet.	MAGLUMI 2000 Plus
4	Pipettor probe	Pipettor	Requires periodic replacement of components and replacement in case of bending damage	1. Using a 1.5 mm wrench, loosen the bolt located directly under the adapter by turning it 2-3 times (it is not necessary to completely unscrew the bolt). Then remove the used probe by pulling it downwards (the probe can be rotated before pulling it out) and then remove the rigid tube connected to the pipettor probe. 2. After removing the used probe, place a new pipettor probe in the rubber sleeve, insert the connecting rigid tube into the new probe, secure it in the adapter and tighten the bolt.	

9.2 Software Maintenance

9.2.1 Software Safety Classification

According to the definitions of software safety classes in IEC 62304 Medical device software - Software life cycle processes, the MAGLUMI® series of automated immunochemiluminescence analyzers is classified as class B (non-serious damage possible) depending on the degree of damage caused by the MAGLUMI. software series.

9.2.2 Software Update Method

Method 1: Reinstall. Uninstall the original software and reinstall the software using a new installation wizard. After the installation is complete, reconfigure the analyzer component parameters.

Method 2: overwrite and update. Close the software, copy and paste the provided patch file into the root directory of the software. Then open the user software. If the software works correctly, it means that the software has been successfully updated.

10 Troubleshooting and Diagnostics

10.1 Common problems

Problem	Possible reasons	Self-Check	Solution
Результаты проверки Abnormal System Test Results: CV and/or RLU out of range	Expired wash liquid, poor water quality or improper preparation of wash liquid	Check wash liquid is too thick, check filth inside the wash tank or check the distilled water conductivity $>1\mu\text{s/cm}$ or not	Clean the wash tank following monthly maintenance, apply distilled water with conductivity $<1\mu\text{s/cm}$; Make system liquid following the instruction.
	Abnormal injection of washer needle	Check the injection volume of 3 pairs of injection needles, see if there is any tubing tangle, air leakage or whether fixing screws of magnetic valve are loose	Contact technical support
	Abnormal soaking Washer hose pump or blocking in washer needles	Check whether the flow speed of the waste tubes is too low or not	Contact technical support
	Contamination inside Chamber or dection window blurred	Check if there are liquid or salt inside chamber	Contact technical support
	The abnormal injection of Starter pump	Check the remaining starters, if there is any strange noise while the starter pumps are working, whether the starter tubes are bent or retracted	Contact technical support
	PMT failure	Check BGW result	Contact technical support
	long using without time maintenance	Check tubing system, if there are dirty or aging	Use Tube cleaning solution to clean tubing system
	Expired starters	Check starter expire date, make sure the storing is under desired condition	Change starters.

Problem	Possible reasons	Self-Check	Solution
High CV of light check	Water quality	Check the distilled water conductivity $>1\mu\text{s/cm}$ or not	Apply distilled water conductivity $<1\mu\text{s/cm}$

Troubleshooting and Diagnostics

	Abnormal cordination of pipettor	Check cordination of pipettors	Contact technical support
	air leakage within pipetting system	Check if there is drop or bubble at the tip of need after washing ; bubbles within pipetting system; leakage or crystal at junction	Check any leakages at junction; Check remaining wash liquid; Check wash position of washing; Re-priming after above checkings, Contact technical support if there is leakage still.
	Sample pump, Washer pump or related valve failure Check warming massage: Initial failed or any liquid leakage or ejection withi	Check warming massage: Initial failed or any liquid leakage or ejection within tubing system	Contact technical support
	Liquid level detection error or suction failure	The results are nearly zero.	Contact technical support
	Dirty damaged pipettor	Check if there are crystals on the surface of pipettor or the Teflon layer damaged	Wipe pipettor if dirty, Contact technical support to replace if damaged
	Contamination inside Chamber or dection window blurred	Check if there are liquid or crystal inside chamber	Contact technical support
	Abnormal injection of Starter pump	Check the remaining starters, if any strange noise whlie working, any bent or retracted happened.	Contact technical support
	PMT failure	Проверьте результаты теста BGW.	Contact technical support
	Long time using without maintenance	Check tubing system, if there are dirty or aging Use Tube cleaning solution to clean tubing sys	Use Tube cleaning solution to clean tubing system.

Problem	Possible reasons	Self-Check	Solution
Test results with high CV	Expired system liquid, poor water quality or improper preparation of wash liquid	Check wash liquid is too thick, check filth inside the wash tank or check the distilled water with conductivity >1μs/cm or not	Clean the wash tank following monthly maintence, apply distilled water with conductivity <1μs/cm; Make system liquid following the instruction of insert

Troubleshooting and Diagnostics

Air leakage within pipetting system	Check if there is drop or bubble at the tip of need after washing; bubbles within pipetting system; leakage or crystal at junction	Check any leakages at junction; Check remaining system liquid; Check wash position of washing; Re-priming after above checkings, Contact technical support if still leakage
Sample pump, Washer pump or related valve failure	Sample pump, Washer pump or related valve failure	Sample pump, Washer pump or related valve failure
Liquid level detection error, suction failure on sample or reagent	Liquid level detection error, suction failure on sample or reagent	Liquid level detection error, suction failure on sample or reagent
Dirty or damaged pipetting needle	Dirty or damaged pipettor	Dirty or damaged pipetting needle
Abnormal injection of washer needle	Check the injection volume of 3 pairs of injection needles, see if there is any tubing tangle, air leakage or whether fixing screws of magnetic valve are loose	Contact technical support.
Abnormal soaking Washer hose pump or blocking in washer needles	Check whether the flow speed of the waste tubes is too low or not	Contact technical support.
Dirty or damaged washer needle or incorrect washer position not	Check the surface of needles: if there are crystals or the Teflon layer damaged Check if the washer lift injection position too low or suction position too high	Wipe pipettor if dirty, Contact technical support
The Chamber or glass window contaminated	Check if have any crystal or liquid in chamber	Contact technical support
Starter injection volume inaccuracy	Check the Starter injection volume, Starter pump, and the Starter tube	Contact technical support
PMT failure	Check the BGW result	Contact technical support
Long time using without maintenance	Check tubing system, if there are dirty or aging.	Use Tube cleaning solution to clean tubing system.
Bubble or others in serum or reagent	Visual inspection.	Remove bubbles
Sample centrifuge insufficient	Visual inspection.	Adjust centrifuge.
Sample tube loaded in rack improper	Visual inspection.	Make the tube touch the bottom of the sample rack.

Problem	Possible reasons	Self-Check	Solution
Cross-Contamination	Probe not clean enough	Check the sample needle position and the contamination outside.	Adjust the sample pipettor position correctly and clean the probe.
	Diluter pump, washer pump or valve abnormal	Check the historical alarm information about the hose pump initialization, check the pump damaged or leakage or not	Contact technical support.
	Probe contaminated or damaged	Check any crystal or Teflon breakage outside of probe	Clean the probe or ask after-sales to replace it
	Wash Concentration expired; storage, the dilute ratio or water improper	Check the Wash Concentration and tank contaminated or not. Check the <i>specific conductance</i> of distilled water less than 1μs/cm or not.	Clean the Wash Liquid Tank, dilute system liquid by pure water(<1μs/cm)
	The volume of the Washer injection abnormal	Check all the Washer injection volume, injection tubes and tighten valves	Contact technical support.
	Hose pump abnormal or waste needle of Washer stuck	Check the flow speed of liquid in waste tube.	Contact technical support.
	Needles of Washer contaminated, damaged or in incorrectly position	Check whether the needles in Washer contaminated or damaged; check the position of washer lift.	Clean the needles in Washer or contact after-sales.
	Chamber or the glass window of chamber lift contaminated Check if have any crystal or liquid in chamber and chamber lift	Check for salt crystals or liquid inside the chamber or on the the chamber lift.	Contact technical support.
	Sample contaminated during the transport or other analyzer	None	Split sample to test or prior at Maglumi system

Problem	Possible reasons	Self-Check	Solution
Low RLU on sample (even 0)	LLD function abnormal; didn't add any sample or reagent	Check the recent results, especially too high or too low RLU sample	Contact technical support.
	Sample insufficient	Visual inspection.	Add sample.

Troubleshooting and Diagnostics

	Starter pump injection abnormally	Check the Starter injection volume, Starter pump and tubes abnormal or not	Contact technical support.
	PMT abnormal	Check the Background in Service.	Contact after-sales.
	Calibration improper, or not validate the calibration curve	Check the Calibration result abnormal or not.	Recalibrate or confirm it.
High RLU on sample	Sample concentration exceed assay range	None	Dilute.
	Wash Concentration expired; storage, the dilute ratio or water improper	Check the Wash Concentration and tank contaminated or not. Check the <i>specific conductance</i> of distilled water less than 1μs/cm or not.	Clean the wash liquid tank, dilute system liquid by pure water(<1μs/cm)
	Washer injection abnormal	Check all the Washer injection volume, injection tubes and tighten valves.	Contact technical support.
	Hose pump abnormal or waste needle of Washer stuck	Check the flow speed of liquid in waste tube.	Contact technical support.
	Needles of Washer contaminated, damaged or in incorrectly position	Check the needles in Washer contaminated or damaged, and the position of washer lift.	Clean the needles or contact after-sales.
	Light leakage of chamber	Check the RLU of other samples in same batch abnormal or not	Contact technical support
	PMT abnormal	Check the Background in Service	Contact technical support
	Blood serum treatment incorrectly	Visual inspection.	Centrifuge the sample in proper condition.

11 Environmental Protection Requirements

The analyzer shall be manufactured and operated in an environmentally friendly manner and in compliance with the rules and requirements of applicable environmental regulations. At the end of its service, the analyzer shall be disposed of in accordance with the instructions governing the disposal of this type of equipment. The analyzer shall be disposed of in accordance with local biosafety regulations. The analyzer may not be disposed along with household waste. The analyzer shall be delivered to a specialized, accredited (licensed) collection point for the recycling of electrical and electronic equipment.

DISPOSAL AND DISCARD

The useful life of the analyzer is 10 years. At the end of its useful life, the analyzer shall be disposed in accordance with the instructions governing the disposal of this type of equipment. If you have any questions regarding disposal, please contact the manufacturer.

Dispose all materials in a safe, acceptable manner that complies with SanPiN 2.1.3684-21 and SanPiN 3.3686-21.

Equipment is disposed by the organizations having a license for this type of activity.

Certain substances in the analyzer shall be disposed subject to the environmental regulations. Disposal should meet the local regulatory requirements.

All departments using analyzers shall dispose liquid waste subject to toxicity regulations and the manufacturer's or authorized dealer's recommendations.

Waste patient liquids with infectious diseases should be disposed using designated disposal devices.

12 List of international standards that apply to the medical product

Item	Standard	Title
1	EN ISO 14971	Medical devices. Application of risk management to medical devices
2	EN 62366-1	Medical devices – Part 1: Application of usability engineering to medical devices
3	IEC/TR 62366-2	Medical devices –Part 2: Guidance on the application of usability engineering to medical devices
4	EN ISO 13485	Medical devices. Quality management systems. Requirements for regulatory purposes
5	EN ISO 18113-1	Diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 1. Terms, definitions and general requirements
6	EN ISO 18113-3	In vitro diagnostic medical devices. Information supplied by the manufacturer (labelling). Part 3. In vitro diagnostic instruments for professional use
7	EN 62304	Medical devices. Software. Life cycle processes
8	EN ISO 15223-1	Medical devices. Symbols to be used with medical device labels, labelling, and information to be supplied. Part 1. General requirements
9	EN 13612	Performance evaluation of in vitro diagnostic medical devices
10	EN 61010-1	Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 1. General requirements
11	EN 61010-2-101	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 2-101. Particular requirements for in vitro diagnostic (IVD) medical equipment
12	EN IEC 61010-2-081	Safety requirements for electrical equipment for measurement, control, and laboratory use. Part 2-081. Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes
13	EN 61326-1	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 1. General requirements
14	EN 61326-2-6	Electrical equipment for measurement, control and laboratory use. EMC requirements. Part 2-6. Particular requirements. In vitro diagnostic (IVD) medical equipment
15	EN IEC 63000	Technical documentation for the assessment of electrical and electronic products with respect to the

13 Storage, Transportation and Packaging

This medical device may be transported by all means of transport in closed type vehicles in accordance with the requirements and rules of transportation of goods applicable to each type of transport.

Disconnect the power cable from the circuit before storage. Protect the analyzer from dust and other contaminants during storage. It is recommended to store the analyzer in its original package. The analyzer should be individually packed before transportation.

Storage/Transportation Requirements

- During transportation, the analyzer shall be kept in a strictly upright position.
- During transportation, the analyzer shall be protected from moisture, water, strong vibration and compression. Handle the analyzer while handling operations.
- The analyzer shall be stored in a well-ventilated area protected from chemicals, corrosive gases, and direct sunlight.
- Storage/transportation temperature -20~55 °C
- Relative humidity during storage/transportation $\leq 93\%$ RH
- Atmospheric pressure during transportation and storage 500 hPa ~ 1060 hPa

14 Manufacturer's Warranty and Liability

The warranty period for the analyzer, its parts and accessories is 12 months

The analyzer manufacturer guarantees its compliance with safety requirements for electronic and medical analyzers. Legislation establishes requirements for installation and operation of the analyzer, which shall be followed by installers and operators.

The manufacturer guarantees the compliance of technical, functional and operational specifications of the analyzers with the requirements of applicable standards and manufacturer's specifications during the lifetime of the analyzer if the user complies with the operating, transportation and storage conditions specified in the operating documentation.

Shenzhen New Industries Biomedical Engineering Co., Ltd.
No. 23, Jinxiu East Highway, Pingshan District, 518122 Shenzhen,
People's Republic of China
Tel: +86-755-21536601 Fax: +86-755-28292740

Authorized representative of the manufacturer in Russia:
OOO CLS (LLC CLS)125363,
Moscow, 6 Novoposelkova St., building 216.1 room 58 Tel: +7-495-9469597.
E-mail: info@cls.msk.com

Appendix A Performance Indicators

A.1 Correctness and reproducibility of sampling

Sample spikes of 5 μL and 200 μL and reagent spikes of 10 μL and 200 μL , respectively, shall conform to the provisions of Table 1.

Table 1 Requirements to correctness and reproducibility for spikes

Nominal spike volume (v)/ μL	Requirement	
	Error	Variation coefficient (VC)/%
5	Max. $\pm 1 \mu\text{L}$	≤ 5
10	Max. $\pm 1 \mu\text{L}$	≤ 5
200	Max. $\pm 1 \mu\text{L}$	≤ 5

A.2 Correctness and fluctuations of the reaction zone temperature control

The temperature offset in the reaction zone shall be within $36.8^\circ\text{C} \pm 0.5^\circ\text{C}$ and fluctuations shall not exceed 0.3°C .

A.3 Light detector area

A.3.1 Analyzer interference

Once the analyzer is turned on, interference with light detection shall not exceed 10%.

A.3.2 Linearity of light values

The linearity of the light values shall be ≥ 0.99 .

A.3.3 Repeatability of luminescence values

When using the luminol method, the coefficient of variation (CV, %) shall be not more than 5%.

A.3.4 Stability of luminescence values

When using the luminol method, the variation of luminol values shall not exceed $\pm 10\%$.

A.4 Contamination transfer rate

Contamination transfer rate shall not exceed 10^{-5} .

A.5 Intra-batch precision of clinical products

Variation coefficient (VC, %) shall comply with Table 2.

Table 2 Requirements to accuracy of clinical articles

Name	Requirements to variation coefficient
HCG (chemiluminescence method)	$VC \leq 5\%$
TSH (chemiluminescence method)	$VC \leq 5\%$

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Instruction for use

Shenzhen New Industries Biomedical Engineering Co., Ltd.

Dear user! Thank you for choosing the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories.

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If you have any questions about the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories, contact your local representative office.

Shenzhen New Industries Biomedical Engineering Co., Ltd.

**Address: No. 23, Jinxiu East Highway, Pingshan District, 518122 Shenzhen,
People's Republic of China(RU)**

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Shanghai International Holding Corp. GmbH (Europe)

Address: Eifestrasse 80, 20537 Hamburg, Germany

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DEVICE INFORMATION

Device Name: Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories

System design: The analyzer consists of a basic setup, accessories, and software. The main installation consists of a material supply module, a liquid module, a temperature control module, a mechanical drive module, an optical signal registration module and an electrical circuit control module.

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Information about this User Manual

Date of publication: 03/2023

Version: 1.0

Required software: ACCORDING to version 1.2.3.0 or other compatible version

Content

WARNINGS	1
DESCRIPTION OF WARNING SYMBOLS	2
PRECAUTIONARY MEASURES	2
PRECAUTIONS FOR USE	5
Chapter 1 MEASUREMENT PRINCIPLE.....	1
1.1 Stages of laboratory analysis	1
1.2 MEASUREMENT PRINCIPLE.....	2
1.3 CALIBRATION METHOD	3
Chapter 2 GENERAL INFORMATION ABOUT THE ANALYZER.....	1
2.1 GENERAL INFORMATION.....	1
2.2 ANALYZER SPECIFICATIONS.....	1
2.3 GENERAL DESIGN OF THE ANALYZER.....	2
2.4 MAIN COMPONENTS OF THE ANALYZER	3
2.4.1 STORAGE COMPARTMENT FOR THE CUVETTE.....	3
2.4.2 THE AREA OF SAMPLES AND REAGENTS	4
2.4.3 THE SAMPLER.....	4
2.4.4 SAMPLE DISK	6
2.4.5 WASHING UNIT.....	7
2.4.6 CHAMBER	8
2.4.7 CUVETTE CAPTURE DEVICE	9
2.5 CONTROL UNIT.....	10
2.6 ANALYZER LABELS.....	10
2.6.1 Module Labels	10
2.6.2 OTHER LABELS.....	15
2.7 REQUIRED MATERIALS (NOT INCLUDED IN THE PACKAGE)	17
Chapter 3 ANALYZER INSTALLATION.....	17
3.1 REQUIREMENTS FOR STORAGE AND TRANSPORTATION OF THE ANALYZER.....	17
3.1.1 Storage requirements	17
3.1.2 Transportation requirements.....	17
3.2 Installation Requirements	17
3.2.1 Requirements for spatial spacing.....	17
3.2.2 Environmental requirements.....	18
3.2.3 Power supply requirements.....	18
3.3 Unpacking and checking.....	18
3.3.1 Unpacking procedure.....	18
3.3.2 Movement and location	19
3.4 Analyzer Installation.....	19
3.4.1 Electrical connection	19
3.4.2 Connecting washing liquid containers.....	19
3.4.3 Connection of liquid waste tanks.....	19
3.4.4 Fuses	20
3.5 Language Setting The default language of the software is English, the user can select the desired	

	language in the user software on the [Settings] page by clicking the <Language> button, then go to the [Language] page and set the desired language.....	20
3.6	Starting and checking the system status	20
Chapter 4	Basic operating procedures.....	21
4.1	Standard operating procedures	21
4.2	Preparation for analysis	22
4.2.1	Pre-launch checks	22
4.2.2	Connecting to the power supply and logging into the software.....	23
4.2.3	Checking the analyzer status	23
4.2.4	Replacement of consumables.....	24
4.2.5	Loading of reagents	26
4.2.6	Uploading samples.....	28
4.2.7	Checking the conditions for performing the analysis	30
4.3	Analyzing	30
4.3.1	Registration of calibration analysis.....	30
4.3.2	Registration of the QC analysis	31
4.3.3	Sample registration	32
4.3.4	Starting the analysis.....	34
4.3.5	Analysis of added samples.....	34
4.4	Monitoring of states during analysis execution	34
4.5	Working with results	34
4.6	Completion of the analysis	35
4.6.1	Power off	35
4.6.2	Actions after power failure	35
4.6.3	Removal of the reagent kit.....	35
4.6.4	Removing and maintaining the sample rack.....	36
4.6.5	Waste disposal	37
Chapter 5	SOFTWARE FUNCTIONS	38
5.1	Description of the software interface.....	38
5.1.1	Main menu.....	38
5.1.2	Interfaces	39
5.1.3	Submenu	39
5.1.4	Toolbar.....	40
5.1.5	Dialog boxes.....	42
5.1.6	Button	42
5.2	Description of [Worklist] screen functions.....	42
5.2.1	Sample registration	42
5.2.2	Test list	44
5.2.3	View details	45
5.3	Description of the screen functions [Result].....	45
5.3.1	Search for results	46
5.3.2	Deleting the results	46
5.3.3	Information about the result.....	46
5.3.4	Repeat measurement.....	47

5.3.5	Sending results.....	48
5.3.6	Recalculation of results.....	48
5.3.7	Exporting results.....	48
5.3.8	Printing the results	48
5.4	Function description of the [Reagent] screen	49
5.4.1	Information about the reagent.....	49
5.4.2	Calibration registration	51
5.4.3	Calibration results.....	51
5.4.4	Information about the reagent.....	52
5.4.5	Reagent Information List.....	54
5.4.6	Reagent consumption statistics.....	55
5.4.7	Off/On reagent.....	56
5.5	Description of the screen functions [QC].....	56
5.5.1	QC results	56
5.5.2	Analysis of control materials	57
5.6	Description of [System Check] screen functions.....	58
5.6.1	System checks.....	58
5.6.2	Loading solution to check the light signal.....	60
5.7	Description of the functions of the [Status] screen.....	60
5.7.1	Consumables.....	61
5.7.2	Environmental parameters (temperature, humidity and voltage).....	64
5.7.3	Wastes.....	64
5.7.4	Analyzer status	66
5.7.5	Shaker status.....	67
5.8	Description of the functions of the [Settings] screen.....	67
5.8.1	Setting Parameters	68
5.8.2	System Settings.....	78
5.8.3	System maintenance	83
5.9	Description of [Home screen] screen functions.....	85
5.9.1	Test states	86
5.9.2	Other status indicators and light indicators.....	86
Chapter 6	SYSTEM MAINTENANCE	87
6.1	Tools.....	87
6.1.1	Accessories (supplied with the analyzer).....	87
6.1.2	User-generated materials	87
6.1.3	Cuvettes	88
6.1.4	Washing liquid.....	88
6.1.5	Starting reagent.....	88
6.1.6	Solution for cleaning system tubes	88
6.2	Daily maintenance	88
6.3	Weekly maintenance.....	89
6.4	Monthly maintenance	89
6.5	Prolonged non-use	90
Chapter 7	The main malfunctions and ways to eliminate them	90

7.1	Classification of faults	90
7.2	Faults not accompanied by an alarm and how to rectify them.....	90
7.3	Faults not accompanied by an alarm and how to rectify them.....	99
Chapter 8	Environmental requirements.....	99
8.1	Disposal and destruction.....	99
Chapter 9	List of international standards applicable to MD.....	100
Chapter 10	Storage, transportation, packaging.....	101
APPENDIX A.	OPERATIONAL CHARACTERISTICS	1
1. 1	Accuracy and repeatability of dosing.....	1
1. 2	Accuracy of temperature control and fluctuations in the reaction area.....	1
1. 3	LC check components.....	1
1. 4	Biomaterial transfer rates.....	1
1. 5	Reproducibility of clinical assays within a single batch	1
APPENDIX B.	NETWORK INTERFACES.....	1
1. 1	Protocol Description.....	1
1. 2	Protocol ASTM	1
1.2.1.	Characteristics of the control code.....	1
1.2.2.	Description of the standard data transfer format.....	1
1.2.3.	Description of the separators	2
1.2.4.	Description of the message types.....	2
1.2.5.	Description of the examples	8
1.3.	HL7 Protocol.....	11
1.3.1.	Description of the standard data transfer format.....	11
1.3.2	Description of the separators	11
1.3.3	Description of the message types	11
1.3.4	Description of message segments.....	13
1.3.5	Description of the examples	19

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This manual provides information on the operation of the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories and describes the design, operation, and use of the analyzer, as well as routine maintenance procedures, troubleshooting, and other aspects of using the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories. The analyzer must be operated in accordance with this manual.

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*Testing: including hormones, tumor antigens, proteins and peptides, liver diseases, cardiomyopathy, myocardial diseases, immune function, autoantibodies, infectious diseases, concentration of vitamins and drugs in the blood, other enzymes, D-dimer, coagulation, microalbuminuria, human immunoglobulin E, procalcitonin.

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DESCRIPTION OF WARNING SYMBOLS

Designation	Text	Description
	Warning	Warnings indicate the presence of a condition in which there is a risk of injury, damage to the analyzer, as well as data errors or biological hazards.
	Note	The note contains important information that you should pay attention to.

PRECAUTIONARY MEASURES

For safe operation of the analyzer, it is necessary to read the manual carefully before starting work. Failure to follow these precautions may result in injury or damage to the analyzer.

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Warning

The user must regularly inspect and maintain the device in strict accordance with the manual in order to avoid disruption of the analyzer and the occurrence of a threat to the safety of personnel.

To ensure safe and reliable operation of the analyzer, installation and maintenance should be performed only by or with the permission of authorized operating engineers and personnel. All analyzer components must be purchased and tested by our company or authorized dealers.

1. WARNING OF INJURIES CAUSED BY MOVING PARTS

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Warning

Do not touch the moving parts of the analyzer and do not interfere with their movement during the operation of the system. Do not place foreign objects in the way of moving parts. Failure to comply with this requirement may result in injury or damage to the analyzer.

To avoid injury or damage to the analyzer, when loading liquid and solid consumables into the analyzer, follow the established rules and properly remove the sealing film, lids or other packaging from the consumables.

2. WARNING OF DANGER OF ELECTRIC SHOCK

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Warning

Make sure that the input voltage meets the requirements for the operation of the analyzer.

Only an authorized service technician may open the front or rear cover of the analyzer while power is applied to the analyzer.

Leakage of reagents, samples and other liquids inside the analyzer may cause malfunction or electric shock. If liquid gets into the analyzer, immediately turn off the power supply of the device and contact the SNIBE maintenance specialists.

If you need to open the front or back cover of the analyzer to replace parts, turn off the analyzer power before opening.

Improper grounding may result in electric shock and damage to the analyzer.

Do not touch items that carry static electricity protection warning labels or perform electrostatic charge removal.

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Warning

Do not use organic solutions for analysis.

The analyzer design is not explosion-proof. Do not use organic solutions with

an ignition temperature of less than 65 °C near the analyzer. Otherwise, it may cause a fire or explosion.

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Warning

Before use, carefully read the safety regulations and precautions given in the safety data sheet for reagents and consumables.

Avoid direct contact of reagents and consumables with hands and clothing. In case of accidental contact of reagents or washing liquid on clothes or skin of hands, wash them immediately with soap and water. In case of accidental contact with the eyes, rinse them immediately with plenty of water and consult an ophthalmologist.

Improper handling of samples can lead to infection. It is forbidden to work with samples, mixtures and liquid waste with bare hands. Contact of these materials with other parts of the body is also not allowed. To avoid infection, always wear gloves, a gauze mask and special clothing while working; if necessary, use protective glasses.

The use of certain types of reagents and washing liquid can lead to skin lesions. When working with reagents and washing liquid, care should be taken to avoid their contact with the skin of hands and clothes. In case of accidental contact of reagents or washing liquid on clothes or skin of hands, wash them immediately with soap and water. In case of accidental contact with the eyes, rinse them immediately with plenty of water and consult an ophthalmologist.

If a small amount of reagent or sample gets on the analyzer during operation, use a clean cotton cloth soaked in alcohol to clean the device before starting work. If a large amount of reagent or sample gets on the analyzer, immediately stop using it and contact an engineer authorized by our company to promptly solve the problem.

Before transporting the analyzer over a long distance, it must be thoroughly disinfected to avoid the risk of infection.

5. PROTECTION FROM LASER RADIATION

Observe the following precautions to prevent burns caused by the barcode scanner's laser radiation.

Warning

Direct exposure of the barcode scanner's laser beam to the retina of the eyes results in eye injury. Avoid direct eye contact with the laser beam of the barcode scanner. Use safety glasses when performing laser beam inspection of a charged barcode scanner.

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Warning

The disposal of certain substances contained in quality control and calibration solutions and liquid waste is subject to contaminant regulations and disposal standards. Observe local waste disposal regulations and requirements specified by the manufacturer or distributor of the respective reagent.

When analyzing samples from patients with infectious diseases, liquid waste should be disposed of by the healthcare facility's infected waste treatment department.

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Warning

Disposal of certain substances contained in spent analyzers is subject to regulations regarding contaminants. Observe the relevant local regulations when disposing of the product.

8. COMPUTER PROTECTION FROM VIRUSES

To protect your computer from viruses, follow the precautions below.

Warning

To prevent damage to the software due to virus infection or misuse, do not use data transfer functions that are not on our authorized list. Sources of spreading computer viruses include floppy disks, U-drives, and networks.

Do not install any other software and hardware other than those specified by our company to avoid incorrect operation of the system software. Do not run other programs while the system is running.

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Note.

Clinical symptoms and test results from other laboratories should be considered when making clinical decisions based on test results.

The User Manual is subject to change without prior notice. If necessary, you

should contact a customer service representative.

The analyzer is intended to be used only by trained personnel who are familiar with basic instrument functions, analysis, and instrument verification.

Do not touch your computer monitor, mouse or keyboard if your hands are contaminated with chemicals.

Do not bend or squeeze the drain tubes. Failure to do so may cause them to become clogged and consequently drain liquid waste from other components, resulting in damage to the analyzer in case of serious irregularities.

The analyzer generates heat during operation, which is expelled through holes in the rear panel of the analyzer. The workroom must be well ventilated to ensure heat dissipation. If necessary, a ventilating device should be used.

However, the analyzer should not be exposed to direct air current. Failure to do so may result in inaccurate analytical results.

Before starting the analyzer for the first time, adjust the settings of all components to ensure the accuracy of their parameters.

There must be no air bubbles in the starting reagents and washing liquid solution in the working concentration. Otherwise, the validity of the analytical results is not guaranteed.

Start reagent containers in the chemiluminescence immunoassay module must be properly connected. Do not use mixed starting reagents. Otherwise, the validity of the analytical results is not guaranteed.

To ensure safe analyzer operation and accurate test results, use compatible cuvettes approved by our company.

To ensure safe analyzer operation and consistent test results, do not use consumables that are past their shelf life.

The analyzer should be turned on at least 30 minutes prior to operation to ensure stabilization of the measurement system.

Before starting the analysis, make sure that the consumables (washing liquid solution at the working concentration, starting reagents, cuvettes, etc.) are sufficient for the analysis.

Before performing any analysis, perform a background wash (BGW) at least once and ensure that the results obtained are within the normal range.

Otherwise, the validity of the analytical results is not guaranteed.

When injecting a liquid sample, make sure that there are no bubbles on the surface of the sample. Failure to do so may result in a sample introduction error. Do not replace or relocate reagents until they are fully incorporated.

Quality control procedures must be followed prior to analyzing the sample. Otherwise, the validity of the analytical results is not guaranteed.

The analyzer is equipped with an alarm indicator light and speaker. If an analyzer error occurs, the alarm light and speaker will sound an appropriate alarm. To ensure proper analyzer operation and accurate test results, take appropriate corrective action immediately after an error occurs. For troubleshooting information, refer to Chapter 7 of this manual: "Basic Troubleshooting and ways to fix them".

Use only SNIBE-approved tubing system washing liquid to clean and maintain the analyzer tubing.

In the event of spillage or ingress of hazardous substances on the surface or in the analyzer housing, properly disinfect the surfaces involved. Do not use cleaners or disinfectants that may react chemically with component materials or parts of the analyzer. If you have any questions regarding the compatibility of a disinfectant or cleaner with the analyzer materials or components, contact the manufacturer or their representative.

In the event of a serious incident involving the use of the device, notify Shenzhen New Industries Biomedical Engineering Co., Ltd. (Snibe) or our authorized representative and the competent authority of the Member State in which you are located.

When you stop using or dispose of the analyzer, the local software data must be deleted.

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3. SYSTEM MAINTENANCE

Warning

Maintenance of the analyzer should be performed on a regular basis following the requirements specified in the appropriate sections of this manual. Improper maintenance can adversely affect the accuracy and reproducibility of analytical results and may result in analyzer malfunction or injury.

The surface of the analyzer may become dusty if the analyzer is not used for a long period of time. To clean the surface, dampen a clean cotton cloth with water, wring it out, and gently wipe the surface of the analyzer. Take the necessary precautions to prevent water from entering the analyzer housing.

Disconnect power to the analyzer before performing maintenance (unless performing live maintenance). Failure to do so may result in analyzer failure or injury.

The analyzer may receive samples from potentially infected patients. Always wear gloves and special clothing during maintenance to prevent infection.

The analyzer does not contain any user-serviceable components. Do not attempt to open the case of the device or remove any components. If necessary, contact an authorized SNIBE service technician for assistance.

6. Contact an authorized service technician to check if the overheat protection system, the liquid level sensor protector, or the entire system is working.
-

4. SAMPLES, REAGENTS AND QUALITY CONTROL SOLUTION

Warning

The presence of substances such as drugs, anticoagulants, and preservatives in the sample may affect the results of the analysis.

The presence of hemolysis, lipemia, and jaundice affects test results.

Make sure there are no clots in the sample. Otherwise, the pipette dispenser will become clogged, seriously affecting the results of the assay.

Samples should be stored properly. Improper storage of samples can lead to changes in their component composition and incorrect analytical results.

To prevent samples from evaporating, do not leave them open for long periods of time. Evaporation of the sample may result in incorrect analytical results.

Improper storage of reagents and quality control solutions can result in incorrect test results and improper analyzer operation, even if these solutions are not expired. Quality control reagents and solutions should be stored and used in accordance with the instructions of the respective reagent manufacturer.

7. After changing the reagent, calibration must be performed. Lack of calibration can adversely affect the accuracy of analytical results.

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6. ELECTROMAGNETIC COMPATIBILITY

Note.

The analyzer complies with the requirements for electromagnetic emission levels and immunity as specified in IEC 61326-2-6.

2. For proper operation of the analyzer, it is the user's responsibility to ensure that the analyzer is operated in accordance with EMC requirements.
3. It is recommended to evaluate the electromagnetic environment before using the analyzer.

Warning

- 1. The analyzer is designed and tested in accordance with the CISPR 11 standard for Class A devices. It is not intended for residential use and may not provide an adequate level of protection against radio emissions in such environments.**

This equipment is intended for use in a professional medical facility. In a home use environment, the device will most likely not work properly. If electromagnetic interference is suspected, proper operation can be restored by increasing the distance between the equipment and the source of the interference.

Do not use this unit close to sources of strong electromagnetic radiation (such as sources of intentional RF radiation), as this may affect the correct operation of the device.

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


Notes

analysis

Fig. 1.1-1. Methodology for performing the analysis

1.2 MEASUREMENT PRINCIPLE

The analyzer, whose operation is based on the principles of immunoreaction, uses heterogeneous magnetic separation technology combined with chemiluminescent reaction to directly label antigens or antibodies with ABEI luminescence labels. This allows to induce a reaction with antigens or antibodies in the analyzed sample so that these antigens or antibodies immobilized by the magnetic particles form an antigen-antibody immune complex. After the free antigen or antibody and other components are removed in the magnetic separation process, a starter reagent is added to the antigen-antibody immune complex to make the tag fluorescent. The luminescence level is then measured from the ROC curve to determine the concentration of the analyte.

The analyzer's photomultiplier tube (PMT) is used to record the light flux emitted during the chemiluminescence reaction. The linear measurement range of the light wavelength recorded by the FEU is 300-650 nm, and the maximum light wavelength emitted from the chemiluminescent reaction is 420 nm. The light signal generated during the chemiluminescence reaction enters the PMT, where it is converted into an electrical signal. Peripheral circuits collect, process, and convert this signal into a digital value representing the magnitude of the light intensity.

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The following describes the process for performing the analysis in the analyzer. After the sample and reagent are added to the cuvette, they undergo three stages of incubation, washing, and separation before the cuvette is placed in the chamber. Start reagent 1 is injected into the cuvette, after 2.5 s start reagent 2 is injected into the same cuvette to start the chemiluminescence reaction. The recording of optical signals by the device begins 0.1 second after the chemiluminescent reaction and lasts for 3.0 seconds, as shown in Figure 1.2-1. The measurement time depends on the incubation time and the order in which the reagent is handled. If the one-step method is used and the incubation time is 1 minute, the measurement time (the period required to obtain the first result) is no more than 7 minutes. If the two-step method is used and the incubation time is 5 minutes, the measurement time (the period required to obtain the first result) will be no more than 11 minutes.

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
Short description:

- The reference curve is calculated based on 10 standard points.
- The values of the relative light units at the two calibration points obtained as a result of calibration are compared with the values of the relative light units of the corresponding concentration on the reference curve.
- The difference between the relative light intensity at the two calibration points and the relative light intensity of the corresponding concentration on the reference curve is then calculated. After that, a linear relationship is established based on the recalculated relative light intensity (Y axis) and concentration (X axis).
- For other points of the reference curve, the difference of the corresponding relative light units is calculated using a compensation curve, after which the values of relative light intensity (Y axis) and concentration (X axis) are recalculated.

The recalculated curve is used as a valid working curve (see Fig. 1.3-1. "The principle of calibration").

.Fig. 1.3-1. The principle of calibration

CHAPTER 2 GENERAL INFORMATION ABOUT THE ANALYZER

2.1 GENERAL INFORMATION

The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories and a set of compatible diagnostic reagents represent an accurate system that allows direct immunochemiluminescent analysis based on the separation of compounds labeled with aminobutyl-ethyl-isoluminol (ABEI) in a magnetic field. This system is used to perform qualitative or quantitative analysis of the content of analytes in samples of human origin. The analyzer automatically performs sample dosing, reagent loading, incubation, rinsing, measurements and calculation of results, which reduces the number of errors in the analysis results, as well as increases their accuracy and reproducibility.

2.2 ANALYZER SPECIFICATIONS

Table. 2.2-1. Analyzer Specifications

Parameter		Specifications
Main characteristics	Test execution speed	200 tests per hour
	Sample type	Serum, plasma, urine, whole blood
	Barcode Type	Code 128, Code 39, Code 93, Codabar, 2/5 Interleaved
	QC	Monthly quality control
Analyzer Characteristics	Sample volume	10-200 µl
	The area of samples and reagents	A total of 72 sample positions and 20 reagent positions on six sample stands
	The sampler	Liquid level detection, liquid level tracking, clot detection, bubble prevention, collision detection, automatic washing
	Reagent volume	10-450 µl
	Temperature in the sample and reagent area	The temperature in the reagent area is 10 °C (±3 °C)
	Number of cuvettes	The total number of loaded cuvettes is 364
	Reaction temperature	37.0 °C ±0.3 °C, temperature fluctuations should not exceed 0.2 °C
	Mixing method	Mechanical vibrations
External interfaces	Interface	Ethernet
Operating conditions	Temperature	10–30 °C
	Relative humidity	≤70 %
	Height above sea level	Not more than 2000 m
	Other	Do not place the device near sources of electromagnetic interference
Storage conditions	Temperature	-20 to 55 °C
	Relative humidity	≤ 93 %
	Atmospheric pressure	50.0–106.0 kPa

The main malfunctions and ways to eliminate them

	Other	Bright sunlight and corrosive gases should not enter the room, It is necessary to ensure good ventilation
Safety category	Degree of protection against electric shock	Type I
	Type of overvoltage	Type II
	Degree of contamination	Degree 2
Analyzer	Weight	132 kg
	External dimensions: length x width x height	90 x 75 x 78 cm
	Package dimensions length x width x height	105.3cm x 88.3 cm x 97.7 cm
	Power supply	100-240 V AC, 50/60 Hz
	Power consumption (V A)	600 V·A
	Noise level	< 70 dB

2.3 GENERAL DESIGN OF THE ANALYZER

The analyzer consists of the main installation, accessories and software. The main installation consists of a material supply module, a liquid module, a temperature control module, a mechanical drive module, an optical signal registration module and an electrical circuit control module.

- The material preparation module consists of a cuvette storage module, a sample and reagent area, a starter reagent module and a container for used cuvettes.
- The liquid module includes modules for the liquid contour of the sampler, the liquid contour of the washing unit, the liquid contour of the chamber, the module of the washing liquid/liquid waste tank, the coolant module of the sample and reagent area and the condensate contour module of the sample and reagent area.
- The temperature control module consists of an incubator heating module (where up to 80 cuvettes can be incubated simultaneously), a sample and reagent area cooling module, and an PMT temperature control module.
- The mechanical drive module consists of a cuvette storage compartment, sample disk, sampler, wash assembly, chamber, gripper and shaker of the sample and reagent area.
- The optical signal registration module consists of a chamber module, a PMT module, and a main control loop module.
- The electrical circuit control module consists of a power supply module, a main control board, a switching panel, a cable system, a number of sensors, motors.
- Software: Automatic immunochemiluminescent analyzer MAGLUMI X3 with accessories, custom version. Release version: V1.


Fig. 2.3-1. Internal design of the analyzer

- | | | |
|--------------------|--------------------------------------|---------------|
| (1) Capture device | (2) Cuvette storage compartment | (3) Chamber |
| (4) Washing unit | (5) Tripod component | (6) Incubator |
| (7) Sample disk | (8) The area of samples and reagents | (9) Sampler |

2.4 MAIN COMPONENTS OF THE ANALYZER

The list of the individual components of each module of the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories is presented below.

2.4.1

STORAGE COMPARTMENT FOR THE CUVETTE

Warning

Do not touch the cuvette storage compartment elevator while the analyzer is operating. Otherwise, injury or damage to the analyzer may occur.

(1) Components


Fig. 2.4-1. Schematic drawing of the cuvette storage compartment

1. Elevator

2. Used unit return device

3. Belt mechanism for
transferring new blocks

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The cuvette storage compartment can be divided into three levels. The lower level stores new blocks moved by a special tape of the transfer mechanism. The middle level stores spent units moved by the return unit. The upper level stores units that are currently in use, moved with a hoist. Two blocks of cuvettes can be stored here.

2.4.2

THE AREA OF SAMPLES AND REAGENTS

(1) Components


Fig. 2.4-2. Schematic drawing of the sample and reagent area

1. Shaker component

2. Door detection sensor

3. RFID tag reader

4. Insulation wool

5. Area component

6. Cover

7. Heat dissipation component

8. Barcode Scanner

(2) Function

Used to house reagent kits and sample racks. A pipette dispenser draws reagents and samples through holes in the panel. The RFID tag reader provides reagent information by reading the information on the tag on the reagent, and the barcode reader provides information on the sample ID and its barcode. The reagent area in the sample and reagent area is cooled with water. The temperature in this area is reduced due to the thermal conductivity properties of the metal from which the component is made, to cool the 20 reagent blocks therein. If it is necessary to shut down the cooling system, disconnect the main power supply.

2.4.3

THE SAMPLER

Warning

Do not touch the sampler while the analyzer is operating. Otherwise, injury or damage to the analyzer may occur.

(1) Components


Fig. 2.4-3. Schematic drawing of the sampler

1. X Axis

2. Y axis

3. Z axis

4. Pipette dispenser

(2) Function

This component is designed to transfer reagents from reagent kits located in the sample and reagent area, and to transfer samples from tubes located on the sample rack to cuvettes. It features liquid level detection, measurement and tracking, clot detection, bubble prevention, collision detection and automatic washing.

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(1) Components


Fig. 2.4-4. Schematic drawing of the sample disk

1. Sample disk agitator

2. Rotating sample disc platform

3. Drive mechanism

(2) Function

This component is designed to move empty cuvettes to the appropriate work areas to load samples and reagents, and to mix components in the cuvettes with added samples or reagents.

2.4.5

WASHING UNIT

**Warning**

Do not touch the washing unit while the analyzer is in operation. Otherwise, injury or damage to the analyzer may occur.


(1) Components

Fig. 2.4-5. Schematic drawing of the washing unit

1. Liquid waste
dispenser

2. Washing liquid
dispensers

3. Washing unit disk
rotating platform

4. Agitator of the
washing unit

(2) Function

The washing unit is primarily designed to remove components of the reaction solution that are not associated with a labeled antigen or antibody of magnetic microparticles in the solid phase by means of three-time washing. The disk of the washing unit moves the cuvettes into the magnetic field and other working areas. Four liquid waste dispensers are designed for aspiration of liquid waste. Three washing liquid dispensers are designed to deliver washing liquid. The agitator of the washing unit is designed to stir the mixture in the cuvettes.

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(1) Components


Fig. 2.4-6. Schematic drawing of the chamber

- | | | |
|--------------------|-------------------------------|------------------------------------|
| 1. PMT module | 2. Tip for starting reagent 2 | 3. Tip for starting reagent 1 |
| 4. Shading modulus | 5. Chamber disc | 6. Liquid waste disposal mechanism |

(2) Function

In the chamber, starting reagents 1 and 2 are injected into a cuvette that has been washed and shaded with a shading module. This exposure of the analyte then results in the generation of an optical signal. This signal is detected, processed and converted into a value representing the light intensity by means of a PMT and a control loop. The liquid waste disposal mechanism is designed to remove liquid waste from cuvettes for which the analysis procedure has been completed.

**Warning**

It is forbidden to touch the cuvette capture device while the analyzer is running. Otherwise, injury or damage to the analyzer may occur.


(1) Components

Fig. 2.4-7. Schematic drawing of the capture device

1. Y axis

2. Z axis

Clamping mechanism

X Axis

(2) Function

The capture device moves the cuvettes between the various components of the analyzer.

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Location:

- 1) One label is on the housing grounding unit.
- 2) One label is on the left support of the reagent area.

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The label is located on the inside panel of the used cuvette container door.

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
The label is located on the inside panel of the starter reagent container door.

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The label is located on the inside panel of the container for starting reagents.

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The label is located on the inside panel of the container for starting reagents.

Laser hazard warning label

Reminds users of the presence of the laser and warns against looking directly at the laser beam.

Located in the laser beam output area of the analyzer: left panel of the reagent area door in the sample and reagent area.

Laser window label

Reminds users where the laser window is located.

Located in the laser beam output area of the analyzer: right panel of the reagent area door in the sample and reagent area.

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Reminds users to connect an external liquid system, high and low concentration liquid waste devices, and indicates the location of the connection.


Located on the left side of the external liquid system connector.

Hint labels for the new and used cuvette block

Remind users of the correct position for loading new cuvettes and returning spent cuvette blocks.


The label is located on the left panel of the cuvette storage compartment door.

Label with information on power switches

Reminds users of the differences between the network interface, main power switch and auxiliary power switch.

Located above the power switch.

On (power)

Off (power)

Cuvette storage label (Cuvette storage compartment)

The cuvette storage compartment is used for storing cuvettes as well as for returning spent cuvette blocks. The label is located on the inside panel of the cuvette storage compartment door.

Cuvette storage

Label on the main power switch on the cover

Reminds users to lock the analyzer during normal operation.

Located under the hinged cover.

Please lock the cover when the analyzer is running. Switch off the main power before opening the cover

Electrical shock warning label

Reminds users not to open the case to avoid electric shock. If maintenance is required, contact our company or an authorized maintenance specialist.

The label is located on the back of the case.

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
Indicates a device designed to record the number of assays.

Located inside, on the right side of the case.

Label warning about the risk of infection

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Location:


- 1) On the front of the used cuvette container.
- 2) On the left side of the external liquid system connector.
- 3) On the inside panel of the sample area door.
- 4) On the cistern for low concentration liquid waste.
- 5) On the cistern for high concentration liquid waste.

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Indicates a cistern that holds high concentration liquid waste.


The label is located on the cistern for high concentration liquid waste.

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Indicates a cistern that holds low concentration liquid waste.

The label is located on the cistern for low concentration liquid waste.

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The label is located on the washing liquid reservoir.

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The label is located on the inside panel of the starter reagent container door.

2.6.2

OTHER LABELS

Label	Description
	Manufacturer
	Authorized Representative in the European Community
	Production date
	Catalog number
	Medical device for in vitro diagnostics
	Serial number
	Note. Refer to the attached documents
	Importer

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
This side up

This symbol indicates that the packaging should be placed vertically during transportation.

Located centered at the top on the packages.

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Located centered at the top on the packages.

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Located centered at the top on the packages.

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
Located centered at the top on the packages.

Stacking height limitation

This symbol indicates that stacking is limited to two packages.

Located centered at the top on the packages.

2.7 REQUIRED MATERIALS (NOT INCLUDED IN THE PACKAGE)

Title	Catalog numbers	NO. RU
In vitro Reagent Kit for automatic analyzers MAGLUMI®	individually	individually
Reaction Cups for automatic analyzers MAGLUMI®	130105000101	/
Wash Concentrate for automatic analyzers MAGLUMI®	130299005M	RZN 2021/14863
System Tubing Cleaning Solution for MAGLUMI® automatic analyzers	130299007M	RZN 2021/14881
Starter reagent kit (Starter 1, Starter 2) for automatic analyzers MAGLUMI®	130299004M	RZN 2021/14858
Luminescence measurement control (Light Check) for MAGLUMI® automatic analyzers	130299006M	RZN 2021/14861
Waste Bag for automatic analyzers MAGLUMI®	130210000201	RZN 2021/14988

Chapter 3 ANALYZER INSTALLATION**3.1 REQUIREMENTS FOR STORAGE AND TRANSPORTATION OF THE ANALYZER****3.1.1 Storage requirements**

The packaged Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories must be stored in well-ventilated rooms in which there is no corrosive gas, at temperatures from -20 to 55 ° C, relative humidity no more than 93% and atmospheric pressure from 50.0 to 106.0 kPa.

3.1.2 Transportation requirements

After packing, the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with accessories should be transported at a temperature of -20 ° to 55 ° and a relative humidity of not more than 93%. Avoid sharp impacts, rain and direct sunlight during transportation.

3.2 Installation Requirements**3.2.1 Requirements for spatial spacing**

To ensure personnel access to the analyzer during operation its installation must be carried out in accordance with the following requirements:

- The distance between the left side of the analyzer and the wall or other devices should be at least 50 cm.
- The distance between the right side of the analyzer and the wall or other devices should be at least 50 cm.
- The distance between the front of the analyzer and the wall or other devices should be at least 125 cm.
- There should be enough space to install a liquid waste drain device.
- There should be enough space near the outlet to properly connect and disconnect the power cord. The analyzer should be installed in such a way that it can be easily disconnected from the power supply.

To ensure maintenance of the analyzer and heat dissipation, observe the following guidelines:

- The distance between the rear panel of the analyzer and the wall or other devices should be at least 50 cm.

3.2.2 Environmental requirements

- Environmental temperature: 10°C - 30°C
- Relative humidity: $\leq 70\%$
- The room should be well ventilated, free of dust, mechanical vibrations, loud noise and interference from the power supply network.
- Do not place the analyzer near an electric brush motor, a source of flashing fluorescent light, and electrical devices that turn on and off frequently.
- Avoid direct sunlight and do not place the equipment near heat sources and air flows. The analyzer should be installed on a flat surface that can withstand a load of at least 1049 kg.

3.2.3

Power supply requirements

- The unit can be connected to an AC power source with the following specifications: 100/240 VAC, 50/60 Hz.
- Rated power of the analyzer: 2200 V·A
- Fuse type: 250 V, 10 A - quick blow fuse.
- Overheating protection device: 2455R-70820258 (Technical specifications: IN THE OPEN STATE — 53 ± 3 ; IN THE CLOSED STATE — $45 \pm 4/-5$, DC 12V, 15A), 9700K51-1638 (technical specifications: IN THE OPEN STATE — 75 °C (± 5 °C), IN THE CLOSED STATE — 50 °C (± 15 °C), 230 V AC current, 13 A).
- For uninterrupted operation of the analyzer, the reserve power of the distribution block must be more than 20%.
- To ensure safety, the analyzer must be connected to a properly grounded power source.
- For safety reasons, an overcurrent protection device with a current rating of 15 A must be installed between the building power supply and the unit's power supply.

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3.3.1

Unpacking procedure

Upon receipt, inspect the packaging of the analyzer and check for external damage. In case of damage, please contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or to a local distributor. After ensuring that the package is intact, open it according to the instructions described below.

- The package must be placed vertically so that the arrow on the label located on the package points upwards.
- Unpack the device and make sure that all components are present according to the packing list. If any component is missing, please contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or to a local distributor.
- Carefully inspect the external surfaces of the analyzer for damage. If damage is detected, contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or your local distributor.

3.3.2**Movement and location**

- For short distance moves, lift the analyzer by directly pulling on the analyzer's rope handles. When using the pen, it is recommended that you perform actions to the right and left of the analyzer. To avoid injury to your hands, do not lift the handle to the edge of the analyzer.
- When moving the analyzer, keep the analyzer upright, do not tilt or place the device on its side. Move the analyzer carefully and slowly.
- Do not allow vibration during handling. All detected faults must be checked and eliminated before moving and operating the device.
- The analyzer should be securely fixed, should not wobble, and the length of the legs should be adjusted so that the analyzer is level.

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3.4.1**Electrical connection**

- Connect the monitor, keyboard, and mouse to the appropriate connectors on the computer.
- Connect the power cables of the system unit and monitor to the appropriate ports.
- Connect the analyzer's power cable to a connector that meets the specified requirements.
- The network interfaces of the analyzer and the computer are connected using a network cable.

3.4.2 Connecting washing liquid containers

- The external washing liquid reservoir, washing liquid reservoir connection port, and washing liquid level sensor connection port are located on the right side of the analyzer front panel.
- The analyzer has two tubes for washing liquid. Insert the connectors of both washing liquid tubes into the appropriate ports on the analyzer; insert the connector of the washing liquid level sensor cable into the appropriate port on the analyzer.
- Then connect the other end of the washing liquid level sensor cable to the pail lid and the other end of the washing liquid tube to the pail socket labeled "System liquid" (Washing liquid).


3.4.3 Connection of liquid waste tanks

- The external liquid waste cistern, liquid waste tank connection port, and liquid waste level sensor connection port are located on the right side of the analyzer's front panel.
- The analyzer has two tubes for liquid waste. Insert the connectors of both liquid waste tubes into the appropriate ports on the analyzer; insert the connector of the liquid waste level sensor cable into the appropriate port on the analyzer.
- Then connect the other end of the liquid level sensor cable to the tank lid and the other end of the liquid waste tube to the tank socket labeled "Thick Waste Liquid" or "Thin Waste Liquid" respectively.

3.4.4 Fuses

Location of 10A fuses,

Fuse specification	Installation location	Location of components
10A	Filter	Analyzer main switch

3.5 Language Setting

The default language of the software is English, the user can select the desired language in the user software on the [Settings] page by clicking the <Language> button, then go to the [Language] page and set the desired language.


Fig.3.5-Language Setting

3.6 Starting and checking the system status

After installing the analyzer, make sure that all power and communication cables between the analyzer

and computer, as well as the washing liquid and liquid waste tubing, are properly and securely connected. Turn on the power and start the system.

1. Turn on the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X3 with Accessories and turn the main and auxiliary switches located on the left side of the analyzer to the on position.
2. After performing a Windows logon, start the user version of the software, log in, and confirm that the parameters of the various analyzer components are correct.
3. After logging into the software and running the user version, the analyzer will enter the initialization state and initialization will be performed for each component.
4. Load consumables such as cuvettes, working concentration washing liquid solution, and starter reagents.
5. Verify that the status of each component displayed in the software is within normal range and that consumables and liquids are sufficient.
6. Completely wash the sampler, washing assembly, and chamber tubes.
7. Follow the system verification procedure. Perform five light signal (LC) checks and ten background wash (BGW) procedures (see Section 4.2.3 Analyzer Condition Check for procedure and expected results). Upon receipt of the system verification results that meet the requirements, the verification is completed and the installation of the analyzer is confirmed.

Chapter 4 Basic operating procedures

This section describes the basic operating procedures of the system. After studying this section, the user will be able to perform standard operating procedures with the system.

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Stages	Procedure
1. Pre-launch checks (1) Checking the power supply (2) Checking the analyzer	Check the power supply for proper operation. Check the appearance and status of the analyzer according to the instructions
2. Connecting the power supply	Power up the analyzer and turn the main and auxiliary switches one after the other to the on position
3. Logging in to the user software	Enter the operator's username and password in the software login dialog box.
4. Checking consumables	Verify that the volume of washing liquid solution at the working concentration, Start Reagent 1, Start Reagent 2, and number of cuvettes are sufficient for the current analysis
5. Checking the analyzer status (1) Status check (2) Confirmation of the results of the system check	On the [Status] screen, confirm the media, consumables, spent media, and shaker status. Wait for all status indicators to turn green before proceeding with the analysis. Follow the system verification procedures and make sure that the measurement values of the LC and BGW tests are within the acceptable range
6. Checking the conditions for	Check the analysis parameters.

The main malfunctions and ways to eliminate them

performing the analysis (1) Checking the analysis parameters (2) Profile settings	Configure the profile according to the current conditions
7. Reagent preparation (1) Checking the remaining reagent volume (2) Checking the mixing time of magnetic microparticles	Make sure that the number of remaining reagent tests is sufficient for the analysis. Make sure that the magnetic microparticles of the reagent have been stirred for at least 30 minutes
8. Calibration of the analysis Registration of the QC analysis	Make sure the reagent is calibrated. Perform calibration for reagents that do not have valid calibration data. Perform the QC analysis
9. Sample registration	Register a standard sample. Register an urgent test (STAT). Register the diluted sample. Register the sample for the QC
10. Starting the analysis	Click the <Start> button to perform the analysis
11. Analysis of the added sample	Register the sample added during the analysis procedure (if the sample is a STAT sample, it must be registered as such; if the sample is a diluted sample, it must be registered as appropriate). Click the <Start> button to perform the analysis
12. Processing the result	Query the analysis results, re-calculate and re-measure, delete or print the results
13. Completion of the procedure Exit from the work program	Log out of the work program
14. Completion of work	Turn off the power to the analyzer and computer
15. After performing the tests	Remove reagent kits and samples from the sample and reagent area, empty the liquid waste tank, used cuvette container, and spent cuvette blocks.

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1. Check the power supply to ensure normal power supply.
2. Check the power and data cables of the analyzer and computer and make sure they are properly and securely connected.
3. Ensure that the washing liquid solution reservoirs at the operating concentration and the starting reagents 1 and 2 are properly connected and contain sufficient liquid.
4. Make sure there are enough cuvettes available.
5. Make sure the liquid waste tank is empty.
6. Ensure that the used cuvette container is empty and the waste cuvette blocks have been emptied.

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4.2.2 software

Connecting to the power supply and logging into the

1. Power up the analyzer and turn the main and auxiliary switches one after the other to the on position
2. After logging into Windows, double-click the user software shortcut on your desktop to launch the software. After starting the software, a login dialog box will appear on the screen. Enter the user name and password and click <Login> to enter the user software interface.
3. After the user logs into the software, the analyzer will enter the initialization state and initialization will be performed for each component.

4.2.3

Checking the analyzer status

4.2.3.1

Status check

On the [Status] screen, check the status of consumables, environmental settings, waste, and shaker. The status of parameters within the normal range is displayed in green.

If the consumables status is displayed as invalid, follow the instructions in the consumables loading information section of this manual.

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If the waste condition is displayed as unacceptable, follow the instructions in the Waste Disposal section of this instruction manual.

If the shaker condition is displayed as unacceptable, contact our service technicians.

**Note.**

After all analyzer components have completed initialization, wait for all status lights to turn green before proceeding with the analysis. If the analyzer does not return to normal operation for an extended period of time, please contact our technical service department.

4.2.3.2

Confirmation of system test results

The system check includes performing a background wash (BGW) and a light check (LC).

Before performing a system check, place the light check (LC) solution in the appropriate positions in the sample and reagent areas. For more information, see 5.6.2 "Loading solution to check the light signal".

Click the [System Check] button on the main menu to go to the [System Check] screen. Click the <Add> button to display the [Registration of data for system verification] screen. Specify the number of BGW

The main malfunctions and ways to eliminate them

checks and LC checks, and then click <Registration> button to register the data for system verification. The [Query] window opens. Check the data for system check, and then click <OK> button to start system check. When finished, check the results of the system check.

Conduct 5 LC checks and 10 BGW checks. The system test is considered passed when results are obtained that meet the following requirements.

BGW	100–1200	coefficient of variation $\leq 10\%$
LC	450000–650000	coefficient of variation $\leq 3\%$

Note.

Analyses may be performed only if the results and coefficients of variation of the BGW and LC checks are satisfactory. It is not recommended to perform sample analysis if the test results are outside the acceptable range (deviating values are highlighted in red in the system test results list). Otherwise, it may adversely affect the validity of the assay results.

The standard LC check should be performed once a week. The LC test should also be performed when analyzer parts are replaced. If no LC test is performed, no liquid load is required for the light test.

BGW washing should be performed every day prior to analysis to ensure accurate results.

4.2.4

Replacement of consumables

If the analyzer does not have enough cuvettes, washing liquid solution at the operating concentration or starting reagents, or if any of the consumables are missing, an alarm will sound. The analyzer will return to normal after the missing consumables have been loaded.

4.2.4.1

Loading cuvettes

- **Loading cuvettes**

New and used cuvette blocks are stored in the cuvette storage compartment. The cuvette storage compartment is located on the left side of the analyzer housing. Open the lower analyzer door and place new cuvette blocks in the new cuvette block transfer area located in the lower compartment to load new cuvette blocks.

The procedure of actions%

1. Log in to the user software. Once the analyzer has completed initialization, check the software for the number of new cuvettes and used cuvettes available.
2. Remove the cuvettes from the packaging and prepare the number of new cuvette blocks required for analysis.
3. Open the bottom door of the analyzer. If there are spent cuvette blocks in the return area, remove them before loading new cuvette blocks.
4. Place the new cuvette assembly in the new cuvette assembly transfer area located in the bottom compartment of the cuvette storage compartment. To complete loading, the new cuvette block will be moved to the new cuvette block use area located in the upper compartment. If the cuvettes at the top of the cuvette block are unused, the block with the newly loaded cuvettes will be temporarily placed at the bottom of the new block

transfer mechanism belt.

Note.

Pay attention to the number of remaining cuvettes before and during analysis.

Ensure that there are no spent cuvette units before loading new units, as cuvette units in the area of use cannot be replaced.

If the analyzer has not been used for two weeks, empty the cuvette storage compartment before use and load a new cuvette assembly to ensure operation with clean cuvettes.

The cuvette blocks and cuvette storage compartment are designed to avoid improper installation. Install these components in the correct direction.

4.2.4.2

Adding washing liquid and starting reagent

- **Adding washing liquid**

The washing liquid reservoir is located on the lower right side of the front panel of the analyzer and is connected to the analyzer with a hose. The remaining amount of the washing liquid can be checked on the [Home Screen] and [Status] screens.

The procedure of actions%

1. Unscrew the cap on the empty washing liquid reservoir that has the liquid level sensor and hose on it. Place the lid in a clean, sealed bag and remove the empty washing liquid container.
2. Remove the pail lid from the sealed bag and place it on the loaded pail of washing liquid. Check that the connection is secure and that there are no liquid leaks.

Note.

The washing liquid solution in the working concentration must be prepared in advance. Ensure that the washing liquid reservoir is free of contaminants.

Pay attention to the volume of remaining washing liquid solution at the working concentration before and during the analysis.

- **Adding the starting reagent**

The procedure of actions%

1. Unscrew the cap of the empty starter container, which has the liquid level sensor and hose on it. Place the lid in a clean sealed bag and remove the empty starter reagent container.
2. Remove the reservoir lid from the sealed bag and place it on a new starter reagent container. Check that the connection is secure and that there are no liquid leaks.
3. Enter 30 in the number of washes in the [Wash] box and click the <Start> button to start the wash.
4. Perform the BGW check using the [System Check] screen.

Note.

Pay attention to the remaining volume of starting reagent before and during the analysis.

**The containers for the starting reagent 1 and 2 as well as the tubes are labeled.
Do not mix up the components during installation!**

4.2.5

Loading of reagents

Scan the RFID tag of the reagent on the left side of the analyzer and place the reagent in the appropriate position in the sample and reagent area. Check the following information on the [Reagent] screen before calibrating, analyzing control materials, and analyzing a sample using a reagent.

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
Fig. 4.2-1. Design of the reagent kit

1. **Sealing film:** serves to seal the reagents. The analyzer sample probe is inserted into the reagent kit through a sealing film to withdraw the reagent.
2. **Handle:** used to load and remove the reagent kit.
3. **RFID tag:** contains information about the reagent kit.
4. **Identification label:** contains information about the components, name, shelf life, etc. of the reagent kit.
5. **Fixing clip:** designed to fix the reagent kit when inserted into the reagent area.
6. **Identification blocker:** designed to lock the button when a reagent kit is inserted into the reagent area for identification.
7. **Mechanism for stirring magnetic microparticles:** together with the shaker of the reagent area, the magnetic microparticles are stirred.

4.2.5.2

Loading the reagent kit

1. **Preparation of the reagent kit**

Ensure that the reagent kit is stored in the correct position (sealing film should be facing up). Do not

invert or shake the reagent kit.

- 1) Remove the reagent kit package from the refrigerator.
- 2) Remove the reagent kit from its packaging.
- 3) Check the sealing film and other elements of the reagent kit for leaks. If there is a leak, contact your distributor as soon as possible.
- 4) Check the liquid in each component of the reagent kit for air bubbles. If air bubbles are present, use a pipette to remove them. Ensure that bubbles are completely absent before use.
- 5) Gently rotate the magnetic microparticle stirring mechanism of the reagent kit to ensure that nothing is obstructing the movement of the mechanism.
- 6) Carefully remove the sealing aluminum film from the reagent kit.
- 7) Remove liquid from the surface of the reagent kit to avoid contamination.

Note.

Remove the sealing aluminum film from the reagent kit before using a new reagent.

Reagents should be prepared, used and stored strictly in accordance with the reagent-specific instructions for use. Do not allow air bubbles to form in the reagents. This may affect the accuracy of reagent collection and analytical results.

Reagents from different kits should not be mixed as this may affect the validity of the assay results.

2. Loading the reagent kit

Before inserting the reagent kit into the reagent area, ensure that the analyzer and software are working properly and that the analyzer has been successfully initialized.

- 1) Open the door of the sample and reagent area.
- 2) Remove unused reagent kits according to the indicator lights of the various channels in the reagent area.
- 3) Holding the reagent kit by the handle, bring the reagent kit with the RFID tag toward the readout area (for approximately 2 seconds). If scanning is successful, a beep will sound.
- 4) Insert the reagent kit along the empty reagent well. When the reagent kit is installed correctly, the analyzer will beep.
- 5) The reagent data will be read by the software. The corresponding item will be highlighted on the [Reagent] screen, where the name of the reagent will also be displayed.
- 6) If an error occurs while reading data, repeat the above steps.
- 7) When the reagent reading is complete, close the sample and reagent area door.

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Indicator light color	Reagent kit status
Doesn't light up	Standby mode, reagent kit can be loaded

The main malfunctions and ways to eliminate them

Lights up green	Standby mode. The reagent kit is installed in a stand
Flashing green	Analysis using the reagent kit is complete, the reagent rack can be removed
Lights up yellow	Startup status. The reagent kit is used
Flashing yellow	Error. For example, the reagent kit was removed during use

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After scanning the RFID tag in the reagent area and properly installing the reagent kit, you should check the number of remaining tests in the reagent information area and verify that the remaining tests are sufficient for this assay. If the reagent volume is insufficient to complete the current sample analysis, install a new reagent kit in the reagent area.

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After the reagent set is installed, the [Reagents] dialog box will display the reagent set name, the number of remaining tests for the reagent, the calibration status, and the magnetic microparticle agitation time in the corresponding reagent rack area. For more information, see 5.4.1 "Information about the reagent".

4.2.6

Loading samples

4.2.6.1

Sample rack design


Fig. 4.2-2. Schematic drawing of a tripod for samples

1. Handle	2. Test tube holder	3. Identification label of the sample stand
4. Sample tube	5. Optocoupler sensor chip	

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- 1) Samples with hemolysis, lipemia, or evidence of microbial contamination should not be used for analysis.
- 2) Ensure that there are no air bubbles in the sample before analysis.

If all requirements are met, place the tube in the sample rack by following the steps below.

- 1) The test tube shall conform to the specified characteristics.
- 2) Carefully insert the test tube into the sample rack.
- 3) If there is a barcode on the tube, it must be turned with the barcode against the open side of the rack so that the barcode scanner can read the required data.

**Note.**

Wear gloves when handling samples to avoid infection.

Make sure the cap of the test tube is open. Failure to do so may result in damage to the pipette dispenser.

If there is insufficient sample volume or no sample, the analysis result still displays the concentration value, but the result contains an "L" label. Thus, if the analyzer gives a warning that the sampler cannot detect the liquid level, another analysis should be performed.

For safety reasons, the sample must meet all listed requirements and performance characteristics to avoid air bubbles and blood clots.

Do not rotate test tubes with the labeled barcode in the sample rack. This may damage the barcode.

After placing the test tube rack in the sample area, do not change the position of the tubes in the rack to avoid confusion of test results.

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- 1) Open the door of the sample and reagent area.
- 2) Remove unused sample racks according to the color of each sample channel indicator in the sample area.
- 3) Select a free sample channel and place the sample rack in it by gently pushing it up so that the barcode on the rack is clearly visible.
- 4) Holding the tripod in a vertical stable position, insert the tripod until the front tripod lock touches the sensor at the end of the sample area and a click is heard.
- 5) If the barcode is not recognized by the system, reload the sample rack to scan the barcode.

The main malfunctions and ways to eliminate them

Indicator light color	Status of the sample rack
Doesn't light up	Standby mode, you can insert a tripod for samples
Lights up green	Waiting
Flashing green	The analyzer has finished loading samples, and the sample rack can be removed
Lights up yellow	The sample rack is loaded
Flashing yellow	Error. For example, the sample rack was removed during use

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4.2.7

Checking the conditions for performing the analysis

Ensure that the parameters and reference range values of the assay are set correctly before performing the assay. Ensure that the assay reagent is calibrated, the calibration curve is correct and is a working curve, and the calibration is valid. Ensure that the profile, dilution, and QC material are set up correctly. (See Section 5.8.1, Parameter Setup, for instructions on setting reagent, profile, dilution, and QC material parameters).

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4.3.1

Registration of calibration analysis

From the main menu, click the <Reagent> button to go to the [Reagent] screen. Select the reagents that require calibration and click the <Calibration registration> button to complete the registration of the calibration assay. Click the <Start> button on the action bar to perform the calibration.


Fig. 4.3-1. Reagent Screen

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4.3.2

Registration of the QC analysis

Click the <Worklist> button on the main menu to go to the [Worklist] screen. Load the sample rack

The main malfunctions and ways to eliminate them

with control materials properly, select the position where the control material is placed in the sample area, and then click the <CC> button to go to the [Data Logging for QC] screen. Select the control material from the list, then select the QC analysis and click <Registration> to register the control material analysis.


Fig. 4.3-3. Registration of the analysis of control materials

4.3.3

Sample registration

Click the <Worklist> button on the main menu to go to the [Worklist] screen.

Load the samples to be analyzed into the sample rack (if there is a barcode on the sample tube, turn the tube with the barcode toward the open side of the rack), then place the rack in the sample area. When a sample rack is placed in the sample area, the scanner reads the sample data and the software will automatically display the read rack ID and barcode.

Select the tripod and position for the samples whose data you want to change, enter the IN of the sample (if you have a barcode, you do not need to enter IN or optional), select the necessary analysis in the analysis selection area, or use the profile function to select several analyses at once and register the samples.

If this sample is to be analyzed urgently (STAT), click on the position of the appropriate STAT sample tube and enter the sample ID (no need to enter the ID if a barcode is present). Then select the STAT assay and click the <STAT> button to prioritize the STAT samples in the assay. If the sample is to be diluted, in the software, click on the position of the desired sample in the rack and then enter the sample ID (no need to enter the ID if a barcode is present). Select the assay to be diluted, click the <Dilute> button and select the dilution factor.

When using the LIS function, you only need to set the settings correctly and connect the LIS server. After entering the sample data, the analysis data for the sample will be automatically received by the analyzer from the LIS server.

Once registration is complete, you can click <Test List> to view the status of all registered tests, delete or change the number of repetitions, and export the data. You can also click the <Overview> button to view the modified analysis data for each sample, the number of registrations for each analysis, and the number of remaining tests for the reagent.

The screenshot displays a software interface for managing a worklist. It is divided into three main sections: 'Sample area', 'Analysis selection', and 'Profile selection'.

Sample area: This section lists six samples (1-6) with IDs NA001 through NA006. Sample 6 is marked with a checkmark. To the right, there is a column for 'C' (Concentration) and a column for 'No.' (Number) with input fields for values 0001 through 0005.

Analysis selection: This section contains a grid of buttons for various tests: Anti-HAV, CMV IgG, CMV IgM, FT3, FT4, HBsAg, HCG/β-HCG, TOXO IgG, TOXO IgM, 25-OH Vit D, AFP, Anti-TPO, CA125, CA153, CA199, CEA, E2, INS, PRL, PROG, and PSA. Navigation arrows are located below the grid.

Profile selection: This section is currently empty, with navigation arrows at the bottom right.

Bottom navigation bar: Includes buttons for 'STAT', 'QC', 'Dilute', 'Test List', and 'Overview'.

Fig. 4.3-4. Worklist

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Note.

Do not open the analyzer cover during operation.

Do not touch or obstruct the pipette, sample disk, or other moving parts while the analyzer is in operation.

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If you need to add analyses during testing, you will need to register the sample before adding it. Select the desired sample rack in the sample area of the [Worklist] screen and the sample for which you want to add an assay, and then register the sample. For more information about the registration methods, see in the section 4.3.3 "Sample registration".

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If you need to add a sample during testing, place the desired sample in the sample rack and properly load the rack into the sample area, then register the sample. For more information about the registration methods, see in the section 4.3.3 "Sample registration".

If there are not enough reagents in the reagent area for additional assays, load the appropriate reagents into the reagent area. If reagents are loaded for additional analyses after registration, click the **<Start>** button on the action bar to complete the continuous loading procedure.

4.4 Monitoring of states during analysis execution

During analysis, the number of samples, tests, consumables, temperature readings, and other analyzer statuses can be monitored in real time from the home screen. If a serious malfunction occurs, click the **<Stop>** button to stop the analysis procedure. For more information about the states, see the 5.7 "Description of the functions of the [Status] screen" section.

4.5 Working with results

Click the **<Result>** button in the main menu to go to the [Result] screen. This screen displays the results of all sample analyses performed during the day.

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To delete results, re-measure, send results over the network, recalculate, and perform other actions on the analysis results, select the checkbox next to the desired result, or simultaneously select multiple results by checking the checkbox in the Sample ID column and clicking the desired button below the results list.

Click the **<Search>** button at the bottom of the [Result] screen to go to the [Result Search] screen. Then specify the search conditions in the criteria field and click **<OK>**. The list of analysis results displays

information about the analysis results from the archive.


Fig. 4.5-1. Results of the analysis

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4.6.1 Power off

When the user has completed all checks and analyses scheduled for the current day, the user can exit the software and disconnect the printer, computer, analyzer, and other devices.

Note.

If the analyzer will not be used for an extended period of time, the main and auxiliary switches should be set to the off position.

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1. Check for unused cuvettes on the horizontal surfaces.
2. Remove all samples from the sample area.
3. Remove all reagents from the reagent area.
4. Empty the liquid waste tank.
5. Empty the used cuvette container and remove the spent cuvette blocks.
6. Check the analyzer panels and ensure that the starter reagent container is free of dirt and stains. If present, remove dirt and stains by wiping surfaces with a clean soft cloth.

4.6.3 Removal of the reagent kit

1. Open the door of the sample and reagent area.

The main malfunctions and ways to eliminate them

2. Holding the handle, remove the reagent kit from the reagent area. If there are no tests left in the reagent kit, dispose of it properly. If reagents are present in the kit, cover the kit cells with additional sealing film to prevent liquid evaporation.
3. Place the reagent kit in the refrigerator, positioning it appropriately.

4.6.4

Removing and maintaining the sample rack

4.6.4.1

Removing the sample rack

1. Open the door of the sample and reagent area.
2. Holding the handle firmly, gently pull the sample rack and remove it completely from the sample area. If there are no samples in the tubes, dispose of the tubes according to the relevant regulatory requirements.
3. If the tubes contain samples that can be used for further analyses, seal them and store them according to the requirements established by the laboratory.

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- **Case 1:** after direct contact with biochemicals.
- **Case 2:** sample tubes cannot be secured in the sample rack.

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1. Take a pointed instrument, pry the tip of the instrument from the inner cell retainer of the sample rack as shown in Figure 4., and gently lift the retainer.
2. Insert the assay tube with the appropriate specifications and check that it is securely fastened.
3. If the tube cannot be secured, follow steps 1 and 2 until the tube is secure.


Fig. 4.6-1. Lifting the retainer

4.6.5

Waste disposal

4.6.5.1

Disposal of used cuvettes

- Disposal of used cuvettes

Used cuvettes are placed in the used cuvette container located in the lower middle section of the analyzer. The software automatically counts the number of spent cuvettes. You can check their status on the [Home Screen] and [Status] screens.

Slide out the cuvette container, remove the waste bag and replace it with a new one, then close the container. Close the waste bag tightly and dispose of it properly.

**Warning**

After emptying the used cuvette container, reset the number of used cuvettes to close the warning message.

Since substances left in the cuvettes can lead to contamination, the cuvettes should be placed in a waste bag. Used cuvettes and liquid waste should be disposed of in accordance with local statutory and regulatory requirements.

Note. The used cuvette container has the following dimensions: 488.5 mm x 143 mm x 145 mm (length x width x height). The user can purchase suitable waste bags himself according to local regulations or with the help of our experts according to clause 2.8.

- Disposal of used cuvette assemblies

Used cuvette assemblies are placed in the waste block return area located in the central compartment of the cuvette storage compartment. See Section 4.2.4.1, "Loading Cuvettes" for details on extraction methods. Used cuvette assemblies do not pose a risk of biological contamination and should be disposed of following standard waste disposal procedures.

4.6.5.2

Liquid waste disposal

The high and low concentration liquid waste outlets are located at the bottom right of the analyzer front panel and are labeled accordingly.

The high concentration liquid waste comes from the sampler and chamber and contains samples, magnetic microparticles, a washing liquid solution at a working concentration, and starting reagents. The concentration of the sample is high. Pollution levels are high.

The low concentration liquid waste comes from the wash unit and contains samples, reagents, and a washing liquid solution at a working concentration. The concentration of the sample is low. Pollution levels are low.

The liquid waste can be stored in a dedicated liquid waste cistern or disposed of in accordance with local legal and regulatory requirements. The waste status can be checked on the [Home Screen] and [Status] screens. If it is indicated that the waste tank is full, waste should be disposed of in a timely manner.

**Warning**

Always wear gloves when disposing of waste.

The liquid waste should be disposed of in accordance with local legal and regulatory requirements.

Refer to Section 6.4 Monthly Maintenance for details on standard procedures for washing the liquid waste tank.

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After logging into Windows, double-click the user software shortcut on your desktop. After starting the software, a login dialog box will appear on the screen. Enter your user name and password, and then click <Log in> to go to the home screen of the user software. The description of the software interface areas is provided below.


Fig. 0-1. Software user interface

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Button	Function
Worklist	Register samples and tests of control materials, run inspections and other procedures.
Result	Check or delete sample results, perform repeat measurements, and other operations.
Reagent	Check reagent information, record reagent calibration, check calibration results, or perform other operations.
QC	View QC results and schedule, and perform other operations.
System check.	View system check information, view system check results, and perform other operations.
Status	Check analyzer temperature and voltage, consumables, etc.
Settings	Configure analysis settings, profile, network, language, etc.
Home Screen	This screen displays the current status of analyzer components, analysis time countdown, etc.

5.1.2

Interfaces

The interfaces (screens) are designed to display the parameters of each of the main functions of the main menu. The name of each individual screen corresponds to the name of the main function in the main menu. The figure below shows the [Worklist] screen.


Fig. 0-3. [Worklist] screen

5.1.3

Submenu

Click the <Menu> button on the main menu to access the submenu of this function. The name of the submenu is the same as the name of the buttons in the main menu. The figure below shows the [Settings] menu.


Fig. 0-4. Menu [Settings]

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-  : <Warning message> button. When an alarm is generated by the analyzer, this button indicates the type of error that occurred. By pressing this button, the user can view the details of the alarm.
-  : <Start> button. When this button is pressed, upon completion of sample registration, calibration, and control material, the software will send a command to the analyzer to start and run tests.
-  - button <Suspend cycle>. If you want to pause the analysis procedure while performing an analysis, click this button and the analyzer will be paused.

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-  - button <Stop>. If you want to stop the analysis procedure while performing the analysis, click this button. The software will then send a stop command to the analyzer and the movement of the components will be stopped immediately.

Note.

Use of this function by the user is not recommended unless special circumstances warrant, such as analyzer malfunction. Once the assay is stopped, all incomplete assay procedures will become invalid.

- - button <Lock screen>. When you click this button the lock screen is displayed, where you can change the user account or lock the software. By entering the user name and password, you can log back into the software. In the locked state, data exchange between the user software and the analyzer will proceed as normal without stopping the analysis procedures. The screen will display the current system time, analysis countdown, analyzer status, number of samples and analyses, and error messages. By tapping the light indicator icon in the lower right corner of the screen, you can mute the alarm speaker.


Fig. 0-5. Lock screen

: button <Exit>. Use this button to exit the software. If this button is pressed during analysis, a message will be displayed indicating that exit cannot be performed at this time.

2. The status area is intended for quick viewing of the current status of some analyzer components.

- TCP/IP status indicator (status of the computer network ports) - the indicator lights green when the network port is successfully enabled and red when the network port is in error.
- PLC status indicator (PLC connection status) - the indicator lights up green when the connection to the PLC is successful and red when the connection fails.
- LIS status indicator (LIS connection status) - the indicator lights up green in case of successful connection to the LIS server and red in case of connection failure.

- Analyzer States - The analyzer can have the following states: Connect, Initiate, Standby, Test, Maintenance, Cycle Suspend, Continue, Stop, Analyzer Failure, and System Check.
- Analysis Countdown - displays the time remaining until the end of the analysis procedure.

5.1.5

Dialog boxes

Dialog boxes open when you click a specific button on the screen. You can open one or more subwindows in the dialog box. The dialog boxes are divided into function windows and information confirmation boxes.

5.1.6

Button

Menus, screens, and dialog boxes contain one or more buttons that are used to launch functions. The three most common buttons are **<OK>**, **<Cancel>** and **<Return>**. The **<OK>** button in the dialog box is used to confirm the operation; the **<Cancel>** button in the dialog box is used to cancel the operation; the **<Return>** button is used to close the dialog box without performing the operation specified in the message of this window.

5.2 Description of [Worklist] screen functions

Click the **<Worklist>** button on the main menu to go to the [Worklist] screen. The [Worklist] function allows you to register assays, check their lists, and perform other actions.

5.2.1

Sample registration











Fig. 5.2-1. Worklist

Decoding of symbols

Area	Symbols	Description
Sample area		The sample rack is not loaded
		The standard sample rack is loaded
		Error loading the sample rack
		Sample rack selected
		The sampling tripod performs sampling
		Standard sample
		STAT Sample
Analysis selection area		The set of reagents for this analysis is installed in the reagent area, the remaining number of tests is greater than zero
		The set of reagents for this analysis is installed in the reagent area, the remaining number of tests is zero
		The reagent kit for this analysis is not installed in the reagent area

5.2.1.1

Registration methods

1. Registration of a single sample

Select the sample rack whose information you want to change, and then enter the sample ID in sequence according to the position of the sample tube in the rack (this step is not required if a barcode is present). Select the sample ID and assay in the assay selection area, or use the profile function to select multiple assays at once and complete the sample registration. If this sample is to be analyzed urgently (STAT), click on the position of the appropriate STAT sample tube and enter the sample ID (no need to enter the ID if a barcode is present). Then select the STAT assay and click the <STAT> button to prioritize the STAT samples in the assay.

If the sample needs to be diluted, select the analysis that requires dilution. Then click the <Dilute> button to select the dilution factor.

The main malfunctions and ways to eliminate them

When using the LIS function, you only need to set the settings correctly and connect the LIS server. After entering the sample data, the analysis data for the sample will be automatically received by the analyzer from the LIS server.

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Select one or more sample racks in the sample area, then select the entire rack and enter the first and last sample ID. After this, click the Enter key to register a series of ID samples on several racks. Then select the analysis or use the profile function to select the analyses to register the samples.

For details on obtaining analysis data using the LIS function, see step 1, "Single Sample Logging".

Sample series logging is not supported for STAT samples.

5.2.2

Test list

Click the <Test List> button on the [Worklist] screen to go to the [Test List] screen. The following functions are available on this screen: display a list of registered sample analyses, change the number of sample repeats, delete already registered sample analyses, and export sample analysis list information in Excel format.

Sample	QC	N°	Tripod position	Sample position	Sample ID	Analysis	Dilution coefficient	Repeats	Status
--------	----	----	-----------------	-----------------	-----------	----------	----------------------	---------	--------

Repeats Delete Export Return

Fig. 5.2-2. Test list

1. Change the number of repeats

You can change the number of repetitions of one analysis using the "+" and "-" buttons in the column indicating the number of repeats. If you want to change the number of repeats for a series of registered assays, select the desired assays, click the <Repeats> button to access the [Change Repeats] screen, and then change the number of repeats using the + or - buttons or manually enter the desired number in the "Change Value" column.


Fig. 5.2-3. Change the number of repeats

2. Deleting the analysis

Select the sample analysis to be deleted and click <Delete>. Only registered analyses can be deleted.

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5.2.3

View details

On the [Worklist] screen, click <Overview> to go to the [Sample Registration Details Overview] screen. In the dialog box that opens, you can view all assays registered for each sample, the number of registrations for each assay, and the number of remaining tests for the reagent.

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Click the <Result> button in the main menu to go to the [Result] screen. The [Result] screen is used to view, search and delete results, view detailed results, repeat measurements, send results over a network, and recalculate, export and print results. This section describes the functions of the [Result] menu sequentially.


Fig. 5.3-1. Results of the analysis

5.3.1

Search for results

On the **[Result]** screen, click the **<Search>** button and set the search criteria on the screen **[Search for a result]** to display the desired results. To display the results of all tests for the current day, simply click the **<Results for this day>** button. The "Records" field in the upper right corner of the screen displays the total number of records in the current list.

5.3.2

Deleting the results

From the list of test results on the **[Result]** screen, select the results you want to delete. Click **<Delete>**, and the **[Query]** dialog box appears. Click the **<OK>** button in the dialog box to delete the entry.

Note.

Please note that deleted test results cannot be recovered.

5.3.3

Information about the result

From the list of test results on the **[Result]** screen, select a result and double-click on it, or click **<Information>** to go to the **[Result Details]** screen and view detailed result information.

Fig. 5.3-2 Information about the result

Click the <Convert to QC> button to go to the [Convert Sample Data to QC Data] screen; in the "Lot No." field, select the lot number of the control material to be converted and click the <Save> button. The results of the analysis will be transferred by the software to the results for the current day for the corresponding batch of control materials.

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Note.

Only samples analyzed on the same day can be re-measured.

Only samples that have been placed in the sample area and not removed can be re-measured.

The remeasurement parameters must match those of the original analysis. If you change analysis parameters, you must register it again.

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5.3.7

Exporting results

In the list of analysis results on the [Result] screen, select the analysis results that you want to export and click the <Export> button to export the results in Excel format.

5.3.8

Printing the results

In the list of test results on the [Results] screen, select the test results that you want to print and click the <Print> button to go to the [Preview]. The user can check the data in the preview mode. After confirmation, click on the print icon in the upper left corner of the dialog box to print the selected test results.


Fig. 5.3-4 Report preview

5.4 Function description of the [Reagent] screen

Click the <Reagent> button in the main menu to go to the [Reagent] screen. The [Reagent] screen functions are used to view reagent and calibration information, log calibration data, and view calibration results and reagent consumption information. This section describes all the functions of the [Reagent] screen in sequence.

5.4.1

Information about the reagent

Using the [Reagent] screen, you can view information about the installation of the reagent kit, as well as the status, name, mixing time of magnetic microparticles and reagent data in the reagent area.


Figure5.4-1. Screen [Reagents]

5.4.1.1

Reagent color marking

Symbols	Description
	The number "1" indicates that the reagent is set to position 1 in the reagent area, and the background color indicates that this reagent has valid calibration data. Each calibration state has its own color marking.  — valid calibration data exists  — there are no valid calibration data  — there are new calibration data  — calibration is in progress  — the calibration period has expired
	"FT3" indicates the name of the current reagent (if a reagent loading error occurs, the message "ERROR" will be displayed here); "99" indicates the number of remaining tests for the current reagent; the icon  indicates that the current reagent is currently in use and the reagent kit cannot be extracted; the background color indicates the current state ("The shelf life of the reagent has not expired, the reagent has not been selected"). Each reagent state has its own color marking:  — the reagent has expired  — the reagent is not installed in the specified position

	 — the shelf life of the reagent has not expired, the reagent has not been selected  — the position is selected
 00:00	"00:00" — countdown of the mixing time of magnetic microparticles of the current reagent
	A check mark indicates that the current reagent is selected

5.4.1.2

Description of the time state

By default, the mixing time of magnetic microparticles is 30 minutes. If the user presses the **<Start>** button on the action bar before the end of the mixing process (before the value "00:00" is displayed), a warning will appear on the screen that the magnetic microparticles are not completely mixed. The analysis can be started only after the user clicks the **<Confirm>** button.

5.4.1.3

The procedure for using reagents

More than one identical set of reagents is loaded into the reagent area. The loaded reagents are used in the following order:

- 1) Reagents are used in ascending order of the batch number.
- 2) Reagents with fewer remaining tests are used.
- 3) Reagents are used from left to right.

5.4.1.4

Reagent data

Basic information about the selected reagent can be viewed in the Reagent Data area.

5.4.2

Calibration registration

On the [Reagent] screen, select the reagents for which calibration is required (multiple reagents can be selected), and click the **<Calibration registration>** button to complete the registration of the calibration analysis. For more information, see 4.3.1 "Registration of calibration analysis".

5.4.3

Calibration results

On the [Reagent] screen, click the **<Calibration Result>** button to go to the [Calibration Result] screen. The calibration results for all reagents will be displayed in the list at the top of the screen. When you select a calibration result, the information corresponding to this result will be displayed in a list at the bottom of the screen. Select the calibration record and click **<Calibr. curve>** to open the [Calibration Log Information] screen and view more detailed information about the calibration result. All the same

The main malfunctions and ways to eliminate them

operations are available for calibration results as for analysis results. For more information, see in the section 5.3 "Description of the functions of the [Result] screen".

No	Reagent	Batch no.	No. of kits	ID of the reference curve	Status	Calibr. 1 — Average	Label	Calibr. 1 — deviation(%)
----	---------	-----------	-------------	---------------------------	--------	---------------------	-------	--------------------------

Calibrator type	No	Status	RLU	Label	Final cycle	End time
-----------------	----	--------	-----	-------	-------------	----------

Buttons: Res. for the day, Search, Calibration curve, Delete, Export, Print, Return

Fig. 5.4-2. Calibration results

5.4.4

Information about the reagent

In the [Reagent] window, double-click the reagent to go to the [Reagent Information] window. The dialog box will display basic information, curve and calibration information, and the corresponding functions for the reagent will be available.

Reagent Name: CMV IgM
ID of the reference curve: 20180716143945
Reagent shelf life: 2022-12-12 00:00:00

No. of the batch: 20191011927
Total tests for reagent: 100
Calibration shelf life: 2022-08-02 16:08:26

No. of kits: 1212
Remaining tests: 96
Date of last calibration: 2022-07-05 16:08:26

Graph: RLU vs Concentration (log scale). The y-axis ranges from 0 to 600,000. The x-axis ranges from 1 to 100. Three curves are shown, all increasing with concentration.

	Calibrator1	Calibrator2
RLU	36052	242726
Average RLU	36052	242726
RLU-CV	0.00	0.00
RLU - deviation(%)	0.0	0.0
Expected concentration	1.530	19.610
Calculated concentration	1.530	19.610
Deviation(%)	0.0	0.0
Permissible concentration deviation	50	50
Permissible deviation of the RLU	50	50

Buttons: Confirm, Reject, Recalculate samples with the same batch no., Recalculate samples with the same kit number, Save, Calibration Log, Print, Return

Fig. 5.4-3 Information about the reagent

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Content	Description
Calibrator 1, Calibrator 2	Results based on two calibration points
RLU	The measured RLU value of the calibrators (up to three measurements for each calibrator)
Average RLU	The average value of the three dimensions of the RLU
RLU – CV	Coefficient of variation for three dimensions of the RLU
RLU – deviation.[%]	The difference (in percent) between the average value of the three RLU measurements and the factory standard value
Expected concentr.	The set calibration concentration value
Calculated concentration	Click the <Calculate Curve> button to calculate the concentration for the last or new calibration curve.
Concentration deviation (%)	The difference (in percentage) between the expected concentration and the actual calibration value
The range of permissible concentration deviations	The permissible range of deviations of the measured calibration concentration values
Permissible deviations of the RLU	The permissible range of deviations of the RLU values during calibration.

5.4.4.1

Operations performed with calibration results

When the calibration analysis is complete, open the [Reagent Information] screen to view calibration and operating curve information. If the calibration result exceeds the acceptable limits, the corresponding calibration parameters will be highlighted in red. The user can evaluate the results and perform the following actions.

1. When an acceptable calibration result is obtained, check the "Confirm" box and click the **<Confirm>** button. The calibration result will be accepted by the system.
2. If you receive an invalid calibration result, click **Reject**, and then click **<Confirm>**. The system will reject and delete the calibration results. After that, the calibration will need to be performed again.
3. If "Recalculate samples with the same batch no." is selected and the **<Confirm>** button is pressed, the latest calibration curve will be used to recalculate the concentration for all results of the current analysis with the same lot number (applicable only to the current day's sample and control material results).
4. If the option "Recalculate samples with the same kit number" has been selected and the **<Confirm>** button is pressed, the latest calibration curve will be used to recalculate the concentration for all results of the current analysis with the same kit number (applicable only to the current day's sample and control material results).
5. To print the current calibration information, click the **<Print>** button.
6. The user can change the values in the "Remaining Tests", "Concentr. deviation tolerance" and "RLU deviation tolerance" fields. The maximum value in the Remaining Tests field is 100.

5.4.4.2

Calibration log

On the [Reagent Information] screen, click the **<Calibration Log>** button to go to the [Calibration Log] screen.

Fig. 5.4-4. Calibration log

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In the calibration log, double-click on an entry or select it and click the **<Details>** button to go to the [Calibration Log Details] screen to view calibration information.

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5.4.5

Reagent Information List

On the [Reagent] screen, click the **<List of Information>** button to go to the [List of Reagent Information] screen. The list of all reagents that have been loaded into the analyzer and a summary of the reagent data is displayed.

To view the desired information, the reagent information can be filtered by specifying the necessary criteria (analysis type, lot number, kit number, or select the "Loaded" option) and then clicking the <Search> button.

Select one or more reagent kit records and click <Delete> to delete the selected records. Click the <Export> button to export the result details in Excel format. Click the <Print> button to print the selected reagent records. Click <Calibration registration> to register the selected reagent kits for calibration.

Fig. 5.4-6. Reagent Information List

5.4.6

Reagent consumption statistics

On the [Reagent] screen, click <Reagent Statistics > to go to the [Reagent Flow Statistics] screen. Statistics on the consumption of reagents used to perform analyses can be viewed for a certain period of time.

Select the start and end time of the period, click <Search> and enter the password of the current account in the window that appears. Then click <OK>, and a list with statistics of reagent consumption for the assays for the specified time period will be displayed.

Select one or more information records and click the <Export> button to export the selected information in Excel format. To print the selected reagent consumption statistics, click the <Print> button.


Fig. 5.4-7. Reagent consumption statistics

5.4.7

Off/On reagent

On the [Reagent] screen, select a reagent kit and click the <Reagent Off/On> button to enable or disable use of the reagent kit.

Select the reagent set to be deactivated and click the <Reagent Off/On> button. The reagent kit will display a "⊖" icon to indicate that it is disabled.

Select a disabled reagent set and click the <Reagent Off/On> button. The shutdown icon will disappear and the reagent kit can be used.

The "Reagent Off/On" function does not support mass deactivation or activation of reagent use.

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5.5.1

QC results

Click on the <QC Result> tab on the [QC] screen to go to the [QC Result] screen. This screen is used to search and delete QC results, view detailed information about QC results, send QC results to the LIS, as well as to recalculate, export and print QC results. All the same operations are available for the results of the system check as for the results of the analysis. For more information, see in the section 5.3 "Description of the functions of the [Result] screen".

Fig. 5.5-1 QC results

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Viewing QC graphs

- QC graphs for the party

At the top of the [QC Analysis] screen, specify the search criteria. The search results will be displayed in the QC results list at the bottom left of the screen. Select the records in the list, after which the corresponding curve for the selected batch will be displayed in the upper area of the QC graph window.

- QC graphs for the month

If it is necessary to check the QC results for one month, when specifying the search criteria it is necessary to check the "Monthly QC" checkbox and select the last or average value as the QC result for the current day. This will allow the QC results for one month to be assessed and a corresponding graph to be plotted. For further steps, see "QC graphs for the batch" above.

From the QC results list at the bottom left of the screen, select the QC message, and then click <QC Data>. The [Ignore QC Data] screen appears with a list of detailed QC results corresponding to the control material lot number. In the list, select the QC results that you want to ignore and click <Ignore>, then the selected results will be recalculated and displayed on the QC graph and in the list with QC conclusions.

Select the QC message from the list of QC results located at the bottom left of the screen and the list of QC findings will be displayed at the bottom right of the screen. The QC conclusions are evaluated in accordance with the QC user rules. For more information, see 5.8.1.10 "Quality control rules".

5.6 Description of [System Check] screen functions

Click the [System Check] button on the main menu to go to the [System Check] screen. The functions on the [System Check] screen are used to perform BGW and LC tests and to view the results of these tests. This section discusses the functions presented on the screen [System Check].

5.6.1

System checks

5.6.1.1

Starting the system check

On the [System Check] screen, click the <Add> button. The [Registration of data for system verification] screen appears as shown below.

On the [Registration of data for system verification] screen, specify the number of BGW and LC checks, and click the <Registration> button to register the system check data.

	Repeats
Washing the sampler	5
Washing of the flushing unit	5
Washing the chamber	5
BGW	10
LC	5

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Note.

During the registration of the system inspection procedure, it is possible to register the inspection procedure for the sampler, the washing unit and the chamber, respectively. However, it is not possible to register the washing procedure itself. In such a case, at least one system test procedure must be enabled to register the washing procedure.

Upon successful registration of the wash procedure, the analyzer will perform this procedure and automatically run a system check upon completion.

After clicking the <Registration> button, the [Query] window will appear with detailed verification data. When the <OK> button is pressed, the analyzer will perform a system check. When the <Cancel> button is pressed, the software will cancel the system check.


Fig. 5.6-2. Query

Note.

Before performing the LC test, place the light signal test solution in the appropriate positions in the reagent area. For more information, see section 5.6.2 "Loading the solution for checking the light signal".

BGW and LC checks must be carried out after the replacement of the starting reagent.

5.6.1.2

System check results

On the screen [System checks] you can search, delete, view information, export, and print system check results. All the same operations are available for the results of the system check as for the results of the analysis. For more information, see in the section 5.3 "Description of the functions of the [Result] screen".

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
Fig. 5.6-3. System check

5.6.2

Loading solution to check the light signal

Light signal test (LC) solution is used to test the accuracy of the analyzer's pipette. During the LC test, the light signal test solution should be set to position 4 of channel 20 of the reagent area.

Strictly follow the instructions for preparing and storing the solution to check the light signal.

5.7 Description of the functions of the [Status] screen

Click the <Status> button on the main menu to access the [Status] screen and related functions (see figure below). This screen displays the status of each analyzer component.

Fig. 5.7-1. Status

The <Status> screen allows you to view the status of analyzer consumables, environmental parameters, and the status of waste and other components.

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If the status of any component corresponds to an alarm value or is abnormal, the <Status> button icon

blinks. To view specific values, click on the corresponding icon to go to the detailed information screen. This section describes all the functions of the [Status] screen sequentially.

5.7.1

Consumables

In the consumables status area of the [Status] screen, you can view status information for cuvettes, washer fluid, starting reagents, and the coolant tank (see figure below). The following provides basic information about the states of consumables.

Fig. 5.7-2 Consumables

5.7.1.1

Cuvettes

Click the <Cuvette> tab on the [Status] screen (see figure below) to access the [Cuvette] screen. This screen is used to display the number of cuvettes in the current cuvette block in real time, as well as the consumption of cuvette blocks in the storage area (in the sliding compartment for storing cuvettes and in the compartment for returning blocks). This screen is also used to indicate the number of cuvettes in the current cuvette block.


Fig. 5.7-3 Cuvettes

1. Description of the cuvettes on the [Status] screen:

When the number of cuvettes is >72, the icon is displayed in green; when the number of cuvettes is between 1 and 72, the icon is displayed in yellow; when the number of cuvettes is 0, the icon is displayed in red.

2. The following storage area states are possible:

Status	Description
	Cuvette blocks are available
	There are no cuvette blocks

3. The following conditions of the compartment for used blocks are possible:

Status	Description
	Cuvette blocks are available
	There are no cuvette blocks

4. The following conditions of the cuvette are possible:

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5. Reset the number of cuvettes.

In the current cuvette area on the right side of the [Cuvette] screen, select the initial cuvette capture position (capture sequence: left to right, top to bottom) and click the button.  The current total number of cuvettes will change accordingly, and the number of cuvettes will reset.

6. Replacing the cuvette block.

In the cuvette block status display area, on the left side of the [Cuvette] screen, click to replace the cuvette block. 

Note.

The number of cuvettes after the reset must correspond to the actual number of cuvettes in the current block.

This function is used only if the cuvette is added manually or if there is no single cuvette. If the initial position of the cuvette capture is selected incorrectly, an error may occur when capturing the cuvette.

5.7.1.2

Washing liquid and starting reagent

Click the "Washing liquid" tab on the [Status] screen to access the [Washing liquid] screen (see figure

below). This screen is used to view the total volume of the washing liquid in real time.

Fig. 5.7-4. Washing liquid

- If the total residual volume is $> 20\%$ of the total possible load volume, the icon is displayed in green.
- If the total residual volume is $< 20\%$ of the total possible load volume, the icon is displayed in green.
- If the total residual volume is $< 20\%$ of the total possible load volume and remains so for a maximum of 20 minutes, the icon is displayed in red.

On the [Status] screen, you can click the [Starting Reagent] tab to access the [Starting Reagent] screen shown in the figure below. This screen is used to view the total volume of starting reagents 1 and 2 in real time.


Fig. 5.7-5. Starting reagent

- If the total residual volume is $> 20\%$ of the total possible load volume, the icon is displayed in green.
- If the total residual volume is $< 20\%$ of the total possible load volume, the icon is displayed in green.
- If the total residual volume is $< 20\%$ of the total possible load volume and remains so for a maximum of 20 minutes, the icon is displayed in red.

5.7.1.3

Cooling tank

In the consumables information area of the [Status] screen, you can view the status of the cooling tank. If the liquid level in the cooling tank is normal, the icon is displayed in green; if the liquid level in the cooling tank is outside the allowable limit, the icon is displayed in red.

5.7.2

voltage)

Environmental parameters (temperature, humidity and

Click on the Environmental Messages area of the [Status] screen to view messages about the temperature and humidity of the analyzer's operating environment, as well as the voltage level of each analyzer component. The following provides basic information about environmental parameters.


Fig. 5.7-6. Environmental parameters

On the [Status] screen, you can click the <T&V> tab to go to the [T&V] screen as shown below. This screen is used to view the temperature or ambient temperature of each component of the analyzer in real time; temperature units can be adjusted. The value displayed on the icon corresponds to the actual value recorded in real time. If the indicator is within the normal range, the icon is displayed in green or yellow; if the indicator is outside the normal range, the icon is displayed in red. Screens with voltage and humidity indicators are arranged in the same way as the screen with temperature indicators.

Fig. 5.7-7. T&V

5.7.3

Wastes

To view the analyzer's liquid and solid waste status messages, click on the waste message area of the [Status] screen (see figure below). The following provides basic information about the waste conditions.


Fig. 5.7-8 Waste

5.7.3.1

Liquid waste

Click on the <Liquid Waste> tab on the [Status] screen to access the [Liquid Waste] screen (see figure below). This screen is used to view the real-time status of the analyzer's liquid waste level (low and high concentration).


Fig. 5.7-9. Liquid waste

If the liquid level is normal, the icon displays in green; if the level is about to reach maximum, the icon displays in yellow; if the level reaches the warning threshold, the icon displays in red.

5.7.3.2

Solid waste

Click the <Solid Waste> button on the [Status] screen (see figure below). Use to view the container states for used cuvettes and empty cuvette blocks in real time.


Fig. 5.7-10. Solid waste

1. If the number of spent cuvettes is less than 400, the icon is displayed in green; if the number of spent cuvettes is between 400 and 520, the icon is displayed in yellow; if the number of spent cuvettes is greater than 520, the icon is displayed in red. You can click the <Reset> button to reset the status of the respective waste materials to zero. The waste bag should then be closed tightly and disposed of properly.
2. When the number of empty cuvette blocks is 0, the icon is displayed in green. When the number of empty cuvette blocks is 1, the icon is displayed in yellow.

Warning

After emptying the waste container, reset the states of the respective waste. Otherwise, the warning message will still be displayed on the screen. Since substances left in the cuvettes can lead to infection, the cuvettes should be placed in a waste bag. Waste should be disposed of in accordance with local legal and regulatory requirements.

5.7.4

Analyzer status

Analyzer disk states can be viewed by touching the Analyzer status information area on the [Status] screen (see figure below). The following presents basic information about the analyzer states.

Fig. 5.7-11 Analyzer states

5.7.4.1

Disk status

Click the <Disk Status> button on the [Status] screen to go to the [Disk Status] screen shown in the figure below. This screen is used to view cuvette positions and real-time status messages for each analyzer disk.

In the central part of the [Disk Status] screen, three disks are displayed (from left to right): the sample disk, the chamber, and the washing unit. The square represents the incubator. The current states of the cuvettes on the disk are indicated by different colors. There is a list of symbols on the left side of the screen. The list of messages for the current analysis performed for the cuvette at the selected position is displayed on the right side of the screen.


Fig. 5.7-12 Disk status

5.7.5

Shaker status

In the shaker status information area of the [Status] screen, you can view the analyzer's shaker status (see figure below). If the speed of the shaker is within the normal range, the icon is displayed in green. If the speed of the shaker is outside the acceptable range, the icon is displayed in red.


Fig. 5.7-13 Shaker status

5.8 Description of the functions of the [Settings] screen

Click the <Settings> button on the main menu to go to the [Settings] screen. The functions of the [Settings] screen are used to configure various parameters and set system values, as well as to perform analyzer maintenance and other actions. This section describes all the basic functions of the [Settings] screen sequentially.


Fig. 5.8-1. Settings

5.8.1

Setting Parameters

The first area shown on the [Settings] screen is the [Parameter Setup] area, which contains the following settings: "Analysis Parameters", "Dilution Setup", "Auto Dilution", "Reflex Testing", "Profile Setup", "Result Setup", "Result Format", "QC Material Producer", "QC Material", and "QC Rules". The following provides basic information about the various parameter settings.

5.8.1.1

Analysis Parameters

Click the <Analysis Parameters.> button to move to the [Analysis Parameters] screen (see figure below). Configure the appropriate parameters for each analysis.

Click the corresponding button in the [Select Analysis] area on the left side of the screen to select the desired analysis. The details of the analysis will be displayed in the [Analysis Details] area.

1. Changing the analysis information

Select the analysis in the [Select Analysis] area, change the required analysis data in the [Analysis Details] area, and click <Save> to complete the details change procedure.

2. Deleting the analysis

In the [Select Analysis] area, select the analysis to be deleted and click <Delete> to complete the deletion procedure.

Setting parameters >> Analysis parameters

Analysis selection

25-OH Vit D	AFP	AI
Alt	ALD	Anti-HAV
Anti-TPO	CA125	CA153
CA199	CEA	CMV IgG
CMV IgM	E2	FT3
FT4	HbAg	HCG/B-HCG
INS	PRL	PROG

Analysis information

Basic information

Name: AI Abbreviation: AI

Channel no.: AI

Unit: ng/ml * 1 = ng/ml

Other settings

☐ Do not automatically load data into LIS if the result is above the normal range

☐ Do not show on other screens

Calibration information

Repeats: 2 Calibration expiration date: 7

Upper limit CV: 15 ☐ General calibration data

Sample information

Repeats: 1

Normal range

0.1 - 6.56

Measuring range

0.1 - 24

Calculation coefficient (Y = aX + b)

a = 1 b = 0

Advanced calculation

☐ Expand Upper limit: 0

Special settings Reference curve Save Delete Return

Fig. 5.8-2. Analysis Parameters

Click the [Reference Curve] button to move to the [Reference Curve] screen (see figure below). The reference curve is the basis of calibration. For each analysis, there is a reference curve that is assigned a specific identification number and batch number.


Fig. 5.8-3. Reference curve

1. Viewing the reference curve

Select an analysis in the [Analysis Selection] area to display the reference curve ID of the corresponding analysis in the [Reference Curve ID] area. Click on the ID of the reference curve from the list to display detailed information about it.

2. Selection of the reference curve

Click on the reference curve in the [ID of the reference curve], change the data, click the <Curve

The main malfunctions and ways to eliminate them

Selection> button to re-fit the reference curve and display the new curve on the screen. Then click the <Save> button to complete the curve selection and save the data.

3. Removing the reference curve

Select the reference curve and click the <Delete> button. Click the <OK> button to delete the reference curve.

4. Printing a reference curve

Select the reference curve and click the <Print> button to print the reference curve and the corresponding analysis data.

Click the <Special Settings> button to access the [Reagent Special Settings] screen shown in the figure below. The <Special Settings> button appears only when AI analysis is selected and allows you to set special settings for the reagent. This button is not shown if other types of analysis are selected.

Special reagent settings

☐ Enable

Interval 500

Range of sample ID for AI, temp. 0 °C

Start 6001

End 6500

Range of sample ID for AI, temp. 37 °C

Start 6501

End 7000

Sample ID range for NA

Start 7001

End 7500

Sample ID range for ALD

Start 7501

End 8000

Information about PRA

Name PRA

Channel no. PRA

Normal range 0.1 ~ 6.56

Abbreviation PRA

Unit ng/mL

Save Cancel

Fig. 5.8-4. Special Settings

To specify whether to use the special settings function, you must select the on/off icon to turn it on or off.

1. Setting the range for the sample ID

Specify the initial value and range of the sample ID for AI, temp. 0 °C. The sample ID ranges for AI, temp. 37 °C, AII and ALD can be set automatically according to the initial value and interval for AI, temp. 0 °C.

2. Information settings for the plasma renin activity (PRA) assay

Basic analysis information on the PRA can be displayed and customized in the analysis information bar. To save the selected settings, click <Save>.

5.8.1.2

Dilution procedure settings

Click the <Breeding Setup> button to go to the [Breeding Setup] screen (see figure below). For analyses that use a sample diluent, up to 10 dilution coefficients can be set. The dilution coefficient can be adjusted by the user.


Fig. 5.8-5. Dilution settings

1. Adding the dilution coefficient

- Select an assay in the [Assay Selection] area, specify a dilution factor in the [Select the dilution coefficient] area to add the corresponding dilution factor to the [Dilution coefficient] area.
- If you want to adjust the dilution factor, select an assay in [Assay Selection], and then click the <Add> button to display the [Change dilution setting] screen. Change the first and second dilution information, then click the <OK> button to add a dilution factor in the [Dilution coefficient] area.


Fig. 5.8-6. Changing dilution settings

2. Changing information about the dilution coefficient

Select an assay in [Assay Selection], select a dilution factor in [Dilution coefficient], and click

The main malfunctions and ways to eliminate them

<Edit> to go to the [Change Dilution Settings] screen to change the settings for the dilution procedure.

3. Removing the dilution coefficient

Select an assay in [Assay Selection] and specify the dilution factor in [Dilution coefficient]. Click the <Delete> button to delete the dilution coefficient.

5.8.1.3

Automatic dilution

To go to the [Auto Dilution] screen (see figure below), click the <Auto Dilution> button. When using the automatic dilution function, a sample whose result exceeds the automatic dilution threshold is automatically registered in the dilution sample area, after which the analyzer dilutes the liquid and completes the analysis.

Setting parameters >> Auto diluting


Fig. 5.8-7. Automatic dilution

1. Adding an automatic dilution procedure

Select an assay in the [Assay Selection] area, and then select a dilution factor in the [Dilution coefficient selection] area (if dilution procedures are not set up for the assay, this area will not contain information; for information on setting up a dilution procedure, see Section 5.8.1.2, "Dilution Procedure Settings"). Enter the threshold value, click the <Add> button, then the new dilution factor will be displayed at the top of the list and the addition procedure will complete.

2. Changing the automatic dilution threshold

Select the assay in the [Assay Selection] area, select the auto dilution factor from the list, enter the critical value, and then click <Save> to complete the procedure for changing the threshold value for auto dilution.

3. Removing the automatic dilution procedure

Select the assay in the [Assay Selection] area, select the auto dilution factor to be deleted from the list, and click <Delete> to complete the deletion process.

5.8.1.4

Reflex testing

Click the <Reflex Test> button to go to the [Reflex Test] screen shown in the figure below. Reflex testing is an analysis that is triggered if the results of the test being performed are in the range that meets the criteria for performing reflex testing.


Fig. 5.8-8. Reflex testing

1. Adding Reflex Testing

Let's say it is necessary for reflex testing B to be performed if the value of the analysis result A is in the range of 50-150. To do this, select assay A in the [Select an assay] area, and select reflex testing B in the [Select an assay to add] area. In the reflex testing settings area B, specify the upper and lower limits (for example, the upper limit is 150 and the lower limit is 50), click the <Add> button, and then add the relevant information to the list of information about reflex testing.

2. Changing the parameters of reflex testing

Select the analysis in the <Select Analysis> area, select the reflex test whose parameters are to be changed in the list of reflex test details, change the upper and lower limits in the corresponding fields, and click <Save> to complete the parameter change process.

3. Removing Reflex Testing

Select an analysis in the <Select Analysis> area, select the reflex test to be deleted in the Reflex Test Details list, and click <Delete> to complete the deletion process.

5.8.1.5

Profile setup

To go to the [Profile Setup] screen (see the picture below), click the <Profile Setup> button. The user can create a profile for combining analyses, which will be displayed in the [Select Profile] area of the worklist. This allows multiple assays to be assigned to a single sample with a single click, simplifying the sample enrollment process.

1. Adding the profile

Select the desired analysis (multiple selections can be made) in the [Select Analysis] area. In the [Information] about profile] field Enter the profile name and click <Add> to complete the addition process.

2. Changing the profile

Select a profile from the [Profile List] list, change the details in the profile details area on the right, and click <Save> to complete the change.

3. Deleting the profile

Select a profile from the [Profile List] list and click the <Delete> button to complete the deletion process.

Setting parameters >> Reflex testing

Analysis selection

25-OH Vit D	AFP	AI
Alt	ALD	Anti-HAV
Anti-TPO	CA125	CA153
CA199	CEA	CMV IgG
CMV IgM	E2	FT3
FT4	HbSAg	HCG/B-HCG
INS	PRL	PROG

Information about reflex testing

NR	Analysis	Lower limit	Upper limit

Lower limit Upper limit

Analysis selection to add

+ Add Save Delete Return

Fig. 5.8-9. Profile setup

5.8.1.6

Format of the results

To go to the [Format of the results] screen (see the picture below), click the [Format of the results] button. The "Format of the results" function is used to determine the number of decimal places in a numeric range.

Specify the upper limit of the numeric range and the corresponding number of decimal places. Specify the upper limit of a numeric range to automatically change the value of the lower limit of the next range. Click the <Save> button to complete the configuration of the result format settings.

Setting parameters >> Result format

Setting the format of the results

Numeric range		Number of digits after comma
0	< 1	3
1	< 10	2
10	< 100	1
100	< Max	0

Save Return

Fig. 5.8-10. Format of the results

Click the <Setting up results> button to go to the [Setting up results] screen shown in the figure below. To customize the information to be displayed in the results, select the appropriate options. The following options are available to customize the results:

Do not display concentration for test results with sampler errors.

Display an out-of-range marker for the result.

Print detailed results when printing sample results, calibrations, QC, and system checks.


Fig. 5.8-11. Setting up results

In the [Result Setup] area, select the desired options and click the <Save> button to complete the modification process.

Click the <Material Manufacturer> button to access the [QC Material Manufacturer] screen shown in the figure below.


Fig. 5.8-12. Manufacturers of QC materials

1. Adding manufacturers of control materials

Enter the manufacturer's name in the Manufacturer Name field and click <Add> to complete the action.

2. Change of information about the manufacturer of control materials

In the list of control material manufacturers, select the manufacturer whose data you want to change,

change the name of the manufacturer in the Manufacturer Name field below and click <Save> to complete the procedure.

3. Removing control material producers

From the list of QC material manufacturers, select the manufacturer you want to delete and click <Delete> to complete the procedure.

5.8.1.9

Control materials

Click the <QC Material> button to move to the [QC Material] screen (see figure below).

Fig. 5.8-13. Materials for QC

1. Operations performed with QC materials

- Click the <Add QC Material> button to open the [QC Material] information area, specify the required material information, and click the <Save> button to complete the addition procedure.
- In the list of control materials, select the control material to be deleted, and then click <Delete QC Material> to complete the procedure of deleting the control material.

2. Operations performed with the analysis of control materials

- From the list of control materials, select the desired control material and click the <Add Analysis> button to open the [Analyze QC] area. Then select the assay to be added to the [QC Analysis] area (multiple assays can be selected). In the [QC Analysis Parameters] area, the newly added assays are displayed in an organized manner and can be switched between using the up and down buttons on the left side of the screen. Fill in the Control Material Analysis Details field as required and click <Save> to complete the procedure for adding the Control Material Analysis.
- In the list of controls, select the control related to the analysis for the control you want to delete, and then in the [QC Analysis] area, select the analysis of the control you want to delete. Click the <Delete Analysis> button to complete the process of deleting the control material analysis.

3. Inventory of control materials

Click the <QC List> button to go to the [QC List] screen. The control material inventory list will display all QC units: [QC material] and [QC analysis].

5.8.1.10

Quality control rules

Click on the <QC Rules> tab to go to the [QC Rules] screen shown in the figure below. The software

The main malfunctions and ways to eliminate them

supports the Westgard rules, allowing the user to set the necessary QC rules for analysis.

Select the analysis in the [Analysis Selection] area, select the QC rules in the [Westgard Rules] area, and click <Save> to complete the QC rules setup.


Fig. 5.8-14. Quality control rules

5.8.2

System Settings

The second area on the [Settings] screen is the [System Settings] area, which contains the following settings: "User Control", "Language", "Network Connection", "Speakers", "Alarm Limit", "General Settings", "Analysis Mode", and "Version Information". The following provides basic information about the various system settings.

5.8.2.1

Version Information

To access the [Version Information] screen (see the figure below), click the <Version Information> button. This screen displays system information and computer version information.


Fig. 5.8-15. Version Information

5.8.2.2

User Management

Click the <User Management> button to display the [User Management] screen shown in the figure below. The software includes user management functions that are used by authorized system administrators

to assign rights to other users. The function of adding new users, changing the name, password and user rights is supported.


Fig. 5.8-16. User Management

1. Adding a user

Enter the information in the [User Information] area, and click <Add> to add a new user.

2. Changing user data

In the [User List] area, select the user whose information you want to change, change the information in the [User Information] area, and click <Save> to complete the user data change procedure.

3. Deleting users

In the [User List] area, select the user whose record you want to delete, and click <Delete> to complete the user deletion procedure.

5.8.2.3

Language

Click the <Language> button to move to the [Language] screen (see figure below). After changing the language, the software will display the text in the appropriate language.


Fig. 5.8-17. Language

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Fig. 5.8-18. Network connection

The software supports two-way data exchange with the LIS server of the medical institution. The analysis information is received from the LIS server of the medical institution (for this purpose it is necessary to manually enter the sample ID or scan the barcode), the analysis results are sent to the LIS server of the medical institution.

1. Network Connection Settings
 - Select the "None" option to operate without using the features available with a network connection.
 - Select the COM port connection and the corresponding COM port number and baud rate in the [COM Port Setting] area to complete the procedure.
 - Select a TCP/IP connection and specify the appropriate IP address and port in the [Network Settings] area.
2. Node setup
 - Analyzer ID: Set the unique channel number of the analyzer.
 - Node ID: configures the unique channel number of the LIS server.
 - Node communication protocol: select the data transfer protocol between the analyzer and the LIS server - "ASTM" or "HL7".

3. Other network connection settings

If necessary, you can set host settings, automatic operations, automatic results download, wait time/repeat, and enable or disable the "Allow counter sum" and "Allow exploration coefficient" options. When all settings are completed, click the <Save> button to complete the network connection settings procedure. After successful connection, the LIS status indicator located on the action bar at the bottom of the screen will turn green.

Click <Speakers> to go to the [Speakers] page (see the figure below). This screen is used to

enable/disable the sound signals through the speakers of the computer or analyzer in case of a warning or error, as well as to select the type of sound signal.


Fig. 5.8-19. Speakers

To go to the [Signal Selection] screen, click the <Selection> button in the [Speaker Setup] area. Then select an audio signal file from the signal list on the [Signal Selection] screen to listen to the signal, and click <Return> to stop listening. Click the <OK> button to complete the beep type selection procedure.

Click the <Save> button on the [Speakers] screen to save the settings.

5.8.2.6

Limitation of the alarm signal

Click the <Alarm Limit> button to go to the [Alarm Limit] screen shown in the figure below. This screen is used to set up alarms that do not need to be set during the analysis.


Fig. 5.8-20. Limitation of the alarm signal

5.8.2.7

General settings

Click the <General Settings> button to go to the [General Settings] screen shown in the figure below. This screen allows you to configure whether a print message dialog box appears when you take a screenshot,

The main malfunctions and ways to eliminate them

as well as enable emergency mode, set the sample tube type, and set the pump speed of the liquid waste washing unit.

No. of sample rack	Tube type
1	
2	
3	
4	
5	
6	

Fig. 5.8-21. General settings

If it is not possible to perform a scan and recognize the installed tripod for samples, you can activate emergency mode. In it, the user can manually edit the sample data to register the analysis.

You can set the type of test tube for the sample and the pump speed of the liquid waste washing unit. The default type of test tube is standard. The pump speed of the liquid waste washing unit is: level 3 > level 2 > level 1.

5.8.2.8 Analysis mode

Click the <Analysis Mode> button to open the [Analysis Mode] screen (see the picture below). It is used to set the priority order of analyses in the system.

Analysis mode
Analysis priority

5.8-22 Analysis mode

1. Analysis mode

Users can select the analysis mode from the drop-down list: "Analysis Priority", "Tripod priority for samples". In this way, you can set up the sequence of sample analyses and QC analyses.

- Analysis priority: Sets of analyses from one batch are sorted according to the analysis category.
 - Sample tripod priority: Test kits from the same batch are sorted according to the sample tripod.
2. Analysis of cross-contamination reagents in priority

If the user has selected the "Tripod first" option, the "Cross contamination reagent analysis priority" option will appear in the interface.

Check this option and all cross-contamination analyses will be sorted in the "Analysis Priority" mode, and the remaining analyses will be performed in the "Sample Tripod Priority" mode.

5.8.3

System maintenance

The third area on the [Settings] screen is the [System Maintenance] area, with which you can perform initialization, washing, daily, weekly and monthly maintenance, as well as troubleshooting components. The following provides basic information about the functions of the system maintenance area.

5.8.3.1

Initialization

To go to the [Initialization] screen, click the <Init.> button. The user can send a command to the analyzer and initialize the analyzer using the screen shown in the figure below.


Fig. 5.8-23. Initialization

1. Sending a command

Enter the desired command in the [Low Level] field, then click the <Send> button to send the command.
2. Initialization

<Init.>: initialization of the analyzer.

<Reset>: Resets the analyzer. When performing this operation, the cuvette is not extracted from the analyzer.
3. Restarting the shaker

Click <Restart Shaker>. This will restart the shaker in the reagent area.

To go to the [Wash] screen (see the picture below), click the <Wash> button.

Maintenance >> Washing

	Repeats	
Washing the sampler	5	☐ Ignore the alarm.
Washing of the flushing unit	5	
Washing the chamber	10	

Start Return

Fig. 5.8-24. Washing

The washing function is used to wash the analyzer tube system. Enter the number of washes of the component and click the <Start> button to finish washing the corresponding tube system of the component. If the analyzer has been installed recently and you need to perform the first wash, click "Ignore. alarm" to avoid stopping washing due to false alarms caused by a large volume of air in the tubes.

Note.

The user must perform the washing before the analyzer starts performing the analysis. In case of replacement of the washing liquid solution in the working concentration, it is necessary to immediately wash the tube system.

5.8.3.3

Daily, weekly and monthly maintenance procedures

Click the <Daily Maintenance> button to move to the [Daily Maintenance] screen (see figure below).

Setting reminders

☒ No reminder

☐ Remind you when the analyzer boots up

☐ Remind you at a set time

06 : 00 : 00

Save

Maintenance procedures

1. Clean the surface of the analyzer with a clean cotton cloth.
- 2.2 Empty the cuvette storage compartments by disposing of the spent consumables they contain.
- 2.1 Empty the used cuvette containers by disposing of the spent consumables they contain.
- 2.3 Empty the high and low concentration liquid waste tanks.
3. Check the quantity of starting reagents. If there are not enough for the next day's tests, refill or replace the reagents.
4. Check the amount of flushing fluid. If there is not enough for the next day's tests, prepare the fluid in advance.
5. Click the button on the right to start the one-button flushing function once.

Start Return

Fig. 5.8-25. Daily maintenance

1. Setting reminders

In the [Setting reminders] area, select the type of reminders. The following options are available: "Do not remind", "Remind when loading", "Remind at the specified time". When selecting the "Remind at a specified time" option, you must specify the time of the daily reminder in the input field. Click the <Save> button to complete the setup.

2. Maintenance procedures

The maintenance procedures must be performed according to the instructions on the screen. If it is necessary to perform a "one-click" cleaning, click the <Start> button for the appropriate procedure, after which the cleaning will be performed.

The weekly and monthly maintenance screens have similar functions. Their detailed description is not provided in this document. The detailed description of the daily, weekly and monthly maintenance procedures is provided in the 0"System Maintenance".

5.8.3.4

Troubleshooting of components

Click the <Troubleshooting Component> button. It must also be pressed in case of an analyzer error to go to the [Troubleshooting] screen (see the picture below).


Fig. 5.8-26. Troubleshooting the component

In the Component initialization area above the dialog box, select one or more components that need to be initialized separately, and then click the <Submit> button to initiate initialization of the selected components.

Click the <Stop> button to stop all tests. Click the <Continue> button to continue performing analyses that can be performed after component initialization is complete.

5.9 Description of [Home screen] screen functions

Click the <Home Screen> button on the main menu to go to the [Home Screen] menu and access its functions as shown in the figure below. This screen is divided into two parts. The central part shows the real-time status of the analysis run, while the right and left parts show the various real-time statuses and lights of the analyzer. The following provides basic information about the components of the home screen.


Fig. 5.9-1. Home Screen

5.9.1

Test states

In the central part of the home screen there is an area displaying the status of analysis execution. It contains the following information (from top to bottom): analyzer status, time remaining for the current analysis, number of samples and analyses.

1. Analyzer states

The Analyzer can have the following states: Connect, Initial, Standby, Test, Maintenance, Cycle Stop, Continue, Stop, Analyzer Failure, and System Check.

2. Remaining time

The software displays the remaining run time of the current analysis in real time. If a new test is loaded during the analysis, the program shows the remaining time until all tests are completed.

3. Samples

The screen shows the sample indicator: number of samples analyzed/total number of samples.

4. Tests

The screen shows the test completion indicator: number of tests completed/total number of tests.

5.9.2

Other status indicators and light indicators

The analyzer component statuses and indicator lights are displayed on both sides of the center area of the home screen. The left side shows the statuses of consumables and waste. If the status of a component is normal, the icon is displayed in green; if the status indicator reaches the notification threshold, the icon is displayed in yellow; if the status indicator reaches the warning threshold, the icon is displayed in red. To view specific states, click on the appropriate icon and a screen with details will be displayed.

The right side of the screen (top to bottom) displays temperature, voltage, calibration pending, time of last BGW check and last maintenance session.

1. Temperature and voltage

If the temperature and voltage values are within the normal range, the corresponding icon is displayed in green; if the temperature and voltage values are outside the normal range, the corresponding icon is displayed in red.

2. Awaiting calibration

If no new calibration is required, the icon is displayed in green; if a new calibration is required, the icon is displayed in red and the reagent analysis for which calibration data are to be confirmed is displayed below the icon.

3. Time of last BGW check

This field displays the time of the last BGW check. If the BGW check was performed on the current day, the icon is displayed in green; if the BGW check was not performed on the current day, the icon is displayed in red.

4. The time of the last service

This field displays the time the last service was performed.

Chapter 6 SYSTEM MAINTENANCE

The operation and scheduled maintenance procedures should be carried out strictly in accordance with the requirements given in the operating instructions to ensure reliable operation, proper operational condition and long service life of the system.

If problems occur during operation that cannot be fixed, as well as problems requiring maintenance and repair that are not covered in this section, contact the maintenance specialists of Shenzhen New Industries Biomedical Engineering Co., Ltd. or to the distributor.

Table 6-1. Information about the components

Code	Name of accessory	Position	Frequency	Note (component replacement)
1	Peristaltic pump tube	Central part of the front right panel of the analyzer	Replace every year or when 73,440 tests are completed and damaged	Contact a service engineer for replacement

6.1 Tools

6.1.1

Accessories (supplied with the analyzer)

Hexagon socket wrench: used to remove the panel and some parts of the analyzer.

Phillips screwdriver: Used to remove and install the tablet and some analyzer components.

Pipette dispenser cleaning tools: used to clean the pipette dispenser.

Positioning tool: Used to calibrate the position of each component.

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Alcohol: for cleaning various parts.

Cotton swab: for cleaning the pipette dispenser.

Liquid waste tank or liquid waste drain tube: used if there are special requirements for draining liquid

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Clamps: for tightening the tubes.

6.1.3 Cuvettes

Cuvettes are used to transfer samples and reagents needed for the reaction. All cuvettes should be purchased only from Shenzhen New Industries Biomedical Engineering Co., Ltd. The use of cuvettes from other manufacturers is prohibited. Our company is not responsible for inaccurate measurements due to the use of cuvettes manufactured by other companies.

6.1.4

Washing liquid

Washing liquid is used to clean all components of the analyzer. Only washing liquids purchased from Shenzhen New Industry Biomedical Engineering Co., Ltd. may be used. If other washing fluids are used, the cuvettes, reagent pipettes, and tubing may not be cleaned sufficiently, and the assay mixture in the cuvette may not be cleaned sufficiently, adversely affecting the accuracy and reproducibility of the assay results. SNIBE assumes no responsibility for inaccurate results when cleaning cuvettes using washing liquid not purchased from SNIBE.

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6.1.6

Solution for cleaning system tubes

The basic composition of the system tube cleaning solution can dissolve protein, which is required to maintain the analyzer in proper condition. Only the system tube cleaning solution purchased from Shenzhen New Industries Biomedical Engineering Co., Ltd. may be used. If another liquid is used to wash the tubing, the effectiveness of the tubing washing may be reduced, which may result in poor accuracy and reproducibility of the assay results. SNIBE assumes no responsibility for inaccurate results when using system tube cleaning solution not purchased from SNIBE.

6.2 Daily maintenance

Materials needed:

Clean cotton cloth

Maintenance Procedures:

1. Use a clean cotton cloth to clean the surface of the analyzer.
2. Empty the used cuvette containers, disposing of the spent consumables they contain.
3. Empty the cuvette storage compartment by disposing of the spent consumables it contains.
4. Empty the high and low concentration liquid waste tanks.

5. Ensure that the volume of starting reagent is sufficient for the next day's assays, add or replace the starting reagent as necessary.
6. Verify that the volume of washing liquid solution at the working concentration is sufficient for the next day's analyses; add washing liquid if necessary.
7. Click the <Start> button on the [Daily Maintenance] screen to have the one-button washing function automatically start once.

6.3 Weekly maintenance

Materials needed:

Cotton swabs, alcohol, clean cotton cloths

Maintenance Procedures:

1. Using an alcohol soaked cotton swab, clean the outside of the pipette dispenser, the outside of the wash needle and liquid waste needle in the wash assembly, and the outside of the chamber probe and liquid waste needle in the chamber.
2. Clean the surface of the cuvette storage compartment with a clean cotton cloth.
3. Click the <Start> button on the [Weekly Maintenance] screen to have the rinse function controlled by one button automatically start twice.

6.4 Monthly maintenance

Materials needed:

Sodium hypochlorite solution with effective chlorine content of 0.05%~0.1% (i.e. 500mg/L~1000mg/L), clean water, system tube cleaning solution, light signal test solution, starting reagent and washing liquid solution in working concentration.

Maintenance Procedures:

1. Replace the starting reagent with 40□ clean water and rinse the chamber 30 times, then replace the clean water with the starting reagent and rinse the chamber another 30 times.
2. Place 1 container of system tubing cleaning solution in the first channel in the reagent area. Use the system tubing cleaning solution to perform the tubing cleaning operation once.
3. Empty the washing liquid tank and the liquid waste tank.
 - 1) Empty the washing liquid tank and the liquid waste tank.
 - 2) Add 2 liters of sodium hypochlorite solution with an effective chlorine content of 0.05%~0.1% (i.e., 500 mg/L~1000 mg/L) to the barrel, cover with a lid, shake back and forth, and leave for 30 min;
 - 3) Pour the cleaned sodium hypochlorite solution and add clean water (2 liters each time) to the cleaning bucket at least 3 times.
4. Fill the bucket with cleaning solution to spray the pipette dispenser and washing assembly, 30 times each;
5. Place the container of light test liquid in the upper right corner of the sample compartment. Perform the BGW test and the light signal test once each. Review the results and record them.
6. System verification requirements and results:

1)	BGW	100–1200	CV ≤ 10 %
----	-----	----------	-----------

2) LC 450000–650000 CV ≤ 3 %

For specific expected values, see the instructions for use of the LC test reagent.

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Replace the washing liquid solution at the operating concentration with clean water.

Click the <Wash> button on the [Settings] screen and wash the sampler, wash assembly, and chamber 30 times each, filling the tubes with clean water.

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1. First, remove the container with the starting reagent and the washing liquid container. Click the <Start> button on the [Wash] screen and wash the sampler, washing assembly, and chamber 30 times each. Then empty the tubes.
2. Add starting reagents 1 and 2.
3. Add a solution of washing liquid at a working concentration.
4. Click the <Start> button on the [Wash] screen and wash the sampler, washing assembly, and chamber 30 times each. Before starting the analysis, make sure that the results of the system check meet the requirements.

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Table 0-1. Faults accompanied by the alarm signal

Faults accompanied by the alarm signal	Analyzer action
Warning alarm	Sample loading and measurements will not be paused, the analyzer will continue to run the assays
Alarm with sample loading stopped	The introduction of the new sample will be stopped, but samples already added to the cuvettes will continue to be analyzed
Alarm with shutdown	All analyzer operations will stop and the analysis results will be unreliable

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Some deviating events may not be accompanied by the alarm. Such events can be detected by checking control procedures, reviewing results, and checking analyzer status.

Table 7.2-1. Faults not accompanied by the alarm and how to rectify them

Signs of malfunction	Possible reason	What needs to be checked	Decision
Abnormal background check result: CV or RLU values outside the permissible range	Either the washing liquid solution at the operating concentration has been stored too long, or a water quality problem in the system or a configuration error has resulted in the washing liquid related malfunction.	Ensure that the washing liquid is not too viscous, that there are no yellow colored contaminants in the washing liquid container, and that the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$	Clean the washing liquid reservoir according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{S}/\text{cm}$, and re-prepare a solution of washing liquid at a working concentration according to the instructions. Rinse the washing unit 30 times.
	Amount of washing liquid supplied to the washing unit is out of the allowable range	Ensure that the three needles of the washing assembly are delivering the normal amount of washing liquid, the supply tube is not bent or leaking, the solenoid valve retaining screws are tight, and the supply tube is not plugged.	Contact technical support
	The washing unit's liquid waste peristaltic pump is malfunctioning or the washing unit's liquid waste removal needle is clogged	Check the liquid flow rate in the liquid waste tube and check for residual liquid in the cuvette.	Contact technical support
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support
	An unacceptable amount of starter reagent was introduced by the starter reagent pump	Ensure that the starting reagent is sufficient, the pump does not produce unusual sounds and the starting reagent tube is not displaced or bent	Contact technical support
	Pollution of the cuvette	Check how long it has been since the cuvettes have been used and make sure they are free of dust and other contaminants	Replace the cuvettes and clean the cuvette storage compartments
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support

The main malfunctions and ways to eliminate them

	The analyzer has been used for a long time and the tubes have not been properly serviced	Check the tubes for signs of wear and dirt	Run the tube washing function once. When performing this procedure, it is necessary to use a washing liquid for the tube system
	The shelf life of the starting reagent has expired	Check the shelf life of the starting reagent and make sure that the storage conditions of the starter reagent meet the necessary requirements	Replace the starting reagent and wash the washing assembly 30 times

Signs of malfunction	Possible reason	What needs to be checked	Decision
High coefficient of variation of LC test results	Problem with water quality in the system	Ensure that the electrical conductivity of the pure water used to prepare the light signal test solution is less than 1 $\mu\text{S}/\text{cm}/\text{cm}$	When preparing the solution, use clean water with an electrical conductivity of less than 1 $\mu\text{S}/\text{cm}/\text{cm}$
	Problem with sampler position	Ensure that the sampler is positioned correctly and is not deflected	Contact technical support
	Air in the pipette pipette tube	Ensure that there are no air bubbles or water droplets on the pipette tip after cleaning, no air bubbles in the sampler tube, and check the tube connections (solenoid valve, piston pump, etc.) for crystallization that may have formed due to leakage	Check the retaining rings at the valve connection for leaks, make sure there is enough washing liquid and that the sampler is in the correct position for washing. Then select the sampler wash option on the Wash screen to perform the wash. If the leak continues, contact technical support
	Sampler piston pump, wash pump, or valve failure	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the washing pump is not damaged and the liquid is not sprayed	Contact technical support
	Liquid level error, LC test solution not added	Check if there are LC test results close to 0	Contact technical support

The main malfunctions and ways to eliminate them

	Pipette dispenser is dirty or damaged	Check the outer surface of the dispenser for crystals or serious damage to the coating	Wipe the pipette dispenser. If there is damage, contact Technical Support for replacement
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support
	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	Pollution of the cuvette	Check how long it has been since the cuvettes have been used and make sure they are free of dust and other contaminants	Replace the cuvettes and clean the cuvette storage compartments
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	The analyzer has been used for a long time and the tubes have not been properly serviced	Check the tubes for signs of wear and dirt	Run the tube washing function once. When performing this procedure, it is necessary to use a washing liquid for the tube system

Signs of malfunction	Possible reason	What needs to be checked	Decision
High coefficient of variation of the analysis results	Either the washing liquid solution at the operating concentration has been stored too long, or a water quality problem in the system or a configuration error has resulted in the washing liquid related malfunction	Ensure that the washing liquid is not too viscous, that there are no yellow colored contaminants in the washing liquid container, and that the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S/cm}$	Clean the washing liquid reservoir according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{S/cm}$, and re-prepare a solution of washing liquid at a working concentration according to the instructions. Rinse the washing unit 30 times

The main malfunctions and ways to eliminate them

	Leakage from the sampler tube	Ensure that there are no air bubbles or water droplets on the pipette tip after cleaning, no air bubbles in the sampler tube, and check the tube connections (solenoid valve, piston pump, etc.) for crystallization that may have formed due to leakage	Check the retaining rings at the valve connection for leaks, make sure there is enough washing liquid and that the sampler is in the correct position for washing. Then select the sampler wash option on the Wash screen to perform the wash. If the leak continues, contact technical support
	Sampler piston pump, wash pump, or valve failure	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the washing pump is not damaged and the liquid is not sprayed	Contact technical support
	Liquid level detection error, sample or reagent not added	Check your most recent test results (as many tests as possible) to see if your the RLU values are too high or too low	Contact technical support
	Pipette dispenser is dirty or damaged	Check the outer surface of the dispenser for crystals or serious damage to the coating	Wipe the pipette dispenser. If there is damage, contact Technical Support for replacement
	Amount of washing liquid supplied to the washing unit is out of the allowable range	Verify that the three needles of the washing assembly are delivering the normal amount of washing liquid, that the supply tube is not bent or leaking, and that the solenoid valve retaining screws are tightened	Contact technical support
	The washing unit's liquid waste peristaltic pump is malfunctioning or the washing unit's liquid waste removal needle is clogged	Ensure that the flexible tube of the liquid waste peristaltic pump in the washing assembly does not need to be replaced; check the liquid flow rate in the liquid waste tube and the presence of residual liquid in the cuvette	Check your test meter, it may be time for a replacement. Contact Technical Support if a replacement or calibration check is required

	Washing unit dispenser dirty or damaged/incorrect cleaning position	Check the outer surface of the dispenser for crystals or other contaminants; make sure the coating on the outer surface has not come off or rusted. Check that the liquid waste removal position and the liquid waste needle position of the washing assembly are correct, and that the needle tube is elastic	Wipe the needle of the washing unit or contact technical support
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support
	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	Pollution of the cuvette	Check how long it has been since the cuvettes have been used and make sure they are free of dust and other contaminants	Replace the cuvettes and clean the cuvette storage compartments
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	The analyzer has been used for a long time and the tubes have not been properly serviced	Check the tubes for signs of wear and dirt	Run the tube washing function once. When performing this procedure, it is necessary to use a washing liquid for the tube system
	Air bubbles or foreign particles were detected in the sample or reagent	External inspection	Remove the air bubbles
	Insufficient centrifugation of serum	External inspection	Standardize the time and speed of serum centrifugation
	The sample tube is not placed in the sample rack or is not in an upright position	External inspection	Install the test tube with the sample for analysis in the sample stand properly

The main malfunctions and ways to eliminate them

Signs of malfunction	Possible reason	What needs to be checked	Decision
Cross contamination	Either the washing liquid solution at the operating concentration has been stored too long, or a water quality problem in the system or a configuration error has resulted in the washing liquid related malfunction	Ensure that the washing liquid is not too viscous, that there are no yellow colored contaminants in the washing liquid container, and that the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$	Clean the washing liquid reservoir according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{S}/\text{cm}$, and re-prepare a solution of washing liquid at a working concentration according to the instructions. Rinse the washing unit 30 times
	Amount of washing liquid delivered to the wash station cuvette is out of the acceptable range	Ensure that the amount of washing fluid supplied to the wash assembly cuvette is normal, that the wash tube is not bent or leaking, and that the solenoid valve retaining screws are tightened	Contact technical support
	The washing unit's liquid waste peristaltic pump is malfunctioning or the washing unit's liquid waste removal needle is clogged	Ensure that the flexible tube of the liquid waste peristaltic pump in the washing assembly does not need to be replaced; check the fluid flow rate of the liquid waste tube	Check the assay counter in the software. If a replacement or calibration check is required, contact your service technician
	Washing unit dispenser dirty or damaged/incorrect cleaning position	Check the outer surface of the dispenser for crystals or other contaminants; make sure the coating on the outer surface has not come off or rusted. Check that the liquid waste removal position and the liquid waste needle position of the washing assembly are correct, and that the needle tube is elastic	Wipe the needle of the washing unit or contact technical support
	Pipette dispenser is dirty or damaged	Check the outer surface of the dispenser for crystals or serious damage to the coating	Wipe the pipette dispenser. If there is damage, contact Technical Support for replacement

	Sampler piston pump, wash pump, or valve failure	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the washing pump is not damaged and the liquid is not sprayed	Contact technical support
--	--	---	---------------------------

Signs of malfunction	Possible reason	What needs to be checked	Decision
The RLU value too low according to the analysis results (or even equal to 0)	Liquid level detection error, sample or reagent not added	Check your most recent test results (as many tests as possible) to see if your the RLU values are too high or too low	Contact technical support
	Insufficient sample quantity	External inspection	Add the sample. If the sample volume is sufficient, contact technical support
	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	BGW value too high	Check to see if the BGW test results for the system prior to analysis were higher than normal	Perform the BGW check again
	Sample concentration is too high	Ensure that the competitive method is applied in the analysis	Perform the dilution
The RLU value too high according to the analysis results	The concentration of the sample exceeds the measurement range	—	Perform the dilution

The main malfunctions and ways to eliminate them

	Either the washing liquid solution at the operating concentration has been stored too long, or a water quality problem in the system or a configuration error has resulted in the washing liquid related malfunction	Ensure that the washing liquid is not too viscous, that there are no yellow colored contaminants in the washing liquid container, and that the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$	Clean the working concentration washing fluid reservoir according to the monthly maintenance procedure, use clean water with an electrical conductivity less than 1 $\mu\text{S}/\text{cm}$, and re-prepare the working concentration washing fluid solution according to instructions
	Amount of washing liquid supplied to the washing unit is out of the allowable range	Ensure that the three needles of the washing assembly supply a constant amount of washing fluid, that the supply tube is not bent or leaking, and that the solenoid valve retaining screws are tightened	Contact technical support
	The washing unit's liquid waste peristaltic pump is malfunctioning or the washing unit's liquid waste removal needle is clogged	Ensure that the flexible tube of the liquid waste peristaltic pump in the washing assembly does not need to be replaced; check the fluid flow rate of the liquid waste tube	Check your test meter, it may be time for a replacement. Contact Technical Support if a replacement or calibration check is required
	Washing unit dispenser dirty or damaged/incorrect cleaning position	Check the outer surface of the dispenser for crystals or other contaminants; make sure the coating on the outer surface has not come off or rusted. Check that the liquid waste removal position and the liquid waste needle position of the washing assembly are correct, and that the needle tube is elastic	Wipe the needle of the washing unit or contact technical support
	Pollution of the cuvette	Check how long it has been since the cuvettes have been used and make sure they are free of dust and other contaminants	Replace the cuvettes and clean the cuvette storage compartments
	Light scattering from the chamber	Check whether the RLU values are abnormally high within a single series of assay results	Contact technical support

	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	Insufficient centrifugation of serum	External inspection	Standardize the time and speed of serum centrifugation

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Point	Standard reference	Title
1	EN ISO 14971	Medical devices Application of risk management to medical devices.
2	EN 62366-1	Medical devices Part 1 Design of medical devices with regard to usability
3	IEC/TR 62366-2	Medical devices Part 2 Guidance on the application of usability design to medical devices
4	EN ISO 13485	Medical devices Quality management systems Requirements for regulatory purposes
5	EN ISO 18113-1	In-vitro diagnostic medical devices, information provided by the manufacturer (labeling), part 1 terms, definitions and general requirements
6	EN ISO 18113-3	In-vitro diagnostic medical devices, information provided by the manufacturer (labeling), part 3 in-vitro diagnostic instruments for professional use
7	EN 62304	Medical devices. Software Life-cycle processes
8	EN ISO 15223-1	Medical Devices Symbols used in the marking of medical devices, on labels and in accompanying documentation Part 1 basic requirements
9	EN 13612	Evaluating the performance of the medical devices for in vitro diagnostics
10	EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1. General requirements
11	EN 61010-2-101	Safety requirements for electrical equipment for measurement, control and laboratory use Part 2-101 Particular requirements for in vitro diagnostic medical equipment
12	EN IEC 61010-2-081	Safety requirements for electrical equipment for measurement, control, and laboratory use-Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes
13	EN 61326-1	Electrical equipment for measurement, control and laboratory use ems requirements part 1: general requirements
14	EN 61326-2-6	Electrical equipment for measurement, control and laboratory use. Electromagnetic compatibility requirements. Part 2-6. Specific requirements for in vitro diagnostic (ivd) medical devices.
15	EN IEC 63000	Technical documentation for the evaluation of electrical and electronic products with regard to the restriction of hazardous substances

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- During transportation, the analyzer must be kept in a strictly upright position.
- During transportation, the analyzer must be protected from moisture, water, strong vibrations and compression. Use caution when loading and unloading the unit.
- The analyzer should be stored in a well-ventilated room protected from chemicals, corrosive gases, and strong sunlight.
- Storage/transportation temperature -20~55 °C
- Relative humidity during storage/transportation $\leq 93\%$ RH

Atmospheric pressure during transportation and storage 500 hPa ~ 1060 hPa

APPENDIX A. OPERATIONAL CHARACTERISTICS

1. 1 Accuracy and repeatability of dosing

When pipetting 10, 200, and 450 μl of sample for analysis (using a single pipette pipette for sample and reagent transfer), the results should meet the requirements presented in Table 1.

Table 1.1-1. Requirements for accuracy and repeatability of dispensing

Nominal dispensing volume (volume)/ μl	Requirements	
	Offset	Coefficients of variation (CV) (%)
10	Within $\pm 1 \mu\text{l}$	≤ 5
200	Within $\pm 5\%$	≤ 2
450	Within $\pm 5\%$	≤ 2

1. 2 Accuracy of temperature control and fluctuations in the reaction area

The temperature in the reaction region should be within the range of $37 \text{ }^{\circ}\text{C} \pm 0.3 \text{ }^{\circ}\text{C}$, with a deviation of no more than $0.2 \text{ }^{\circ}\text{C}$.

1. 3 LC check components

1.3.1. Analyzer noise level

During operation of the analyzer when the stable operating condition is reached, the noise level should be in the range of 100-1200, the coefficient of variation should be no more than 10%.

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1. 4 Biomaterial transfer rates

The transfer index should be $\leq 10^{-5}$.

1. 5 Reproducibility of clinical assays within a single batch

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APPENDIX B. NETWORK INTERFACES

1. 1 Protocol Description

The software utilizes ASTM E1394 and HL7 data transfer protocols.

All data transferred between the software and the LIS system is stored in the logs/date/online folder.

1. 2 Protocol ASTM

1.2.1. Characteristics of the control code

Table 1.2-1. Characteristics of the control code

Code	Description	Corresponding ASCII code
<ENQ>	Query	0x05
<ACK>	Confirmation of response	0x06
<STX>	Start of the document	0x02
<ETX>	End of the document	0x03
<CR>	Carriage return	0x0D
<EOT>	Completing data transfer	0x04

1.2.2. Description of the standard data transfer format

<ENQ><STX>content<STX><EOT>

For example:

<ENQ><STX>TEST<STX><EOT>

Table 1.2-2. Description of the standard data transfer format

Code	Description	Corresponding ASCII code
<ENQ>	Query	0x05
<STX>	Start of the document	0x02
T	The letter T	0x54
E	The letter E	0x45
S	The letter S	0x53
T	The letter T	0x54
<ETX>	End of the document	0x03
<EOT>	Completing data transfer	0x04

1.2.3. Description of the separators

The ASTM E1394 protocol provides 4 transmission data separators. For more information, see section A.4 of the E1394 protocol. The value of each separator is shown in the table below.

Table 1.2-3. Description of the separators

Separators	Description	Corresponding ASCII code
<CR>	Record separator	0x0D
	Field separator	0x7C
\	Repeat separator	0x5C
^	Component separator	0x5E
&	Special separator	0x26

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1.2.4.1 Recording the message header (H)

Note.

The information in this section is provided in accordance with Section 7 of ASTM Protocol E1394.

The message header entries must be at the beginning of all transmission entries. They contain basic information about protocols and basic data transmission.

Example:

H|^&||PSWD|MAGLUMI X3|||||LIS||P|E1394-97|20180717<CR>

Table 1.2-5. Recording the message header

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	7.1.1	Message type	H	1	Yes
2	7.1.2	Separator definition	\ ^ &	4	Yes
4	7.1.4	Password	PSWD	20	No
5	7.1.5	Sender name	MAGLUMI X3	20	Yes
10	7.1.10	Receiver name	LIS	20	Yes
12	7.1.14	Processing method	P	1	Yes
13	7.1.13	Protocol version number	E1394-97	10	Yes
14	7.1.14	Date	YYYYMMDD	14	Yes

Note.

The total number of fields in the protocol is 14. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software. The "|" separator should also be used to separate different fields if some of these fields do not require values.

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1.2.4.2 Recording of the patient information (P)

Note.

The information in this section is provided in accordance with Section 8 of ASTM Protocol E1394. Information related to the patient should be presented using this type of record.

Example:

P|1|||ABC|||F<CR>

Table 1.2-6. Recording patient information

Field	E1394	Field name according to ASTM	Content	Maximum length	Necessarily
-------	-------	------------------------------	---------	----------------	-------------

		protocol			
1	8.1.1	Message type	P	1	Yes
2	8.1.2	Serial number	1	6	Yes
6	8.1.6	Patient Name	ABC	30	No
9	8.1.9	Sex	M, F, U	1	No

Note.

Since only the first and second fields are mandatory, the patient information record is usually presented in the format "P|1 <CR>".

The total number of fields in the protocol is 35. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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1.2.4.3 Recording analysis information (O)

Note.

The information in this section is provided in accordance with Section 9 of ASTM Protocol E1394.

The information for each analysis should be presented using a record of this type.

Example:

O|1|1234567||^^ALT^^1:50|R<CR>

Table 1.2-7. Recording analysis information

Table 12-7 Recording analysis information						
Field		E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1		9.1.1	Message type	O	1	Yes
2		9.1.2	Serial number	1	6	Yes
3		9.1.3	Sample ID	1234567	22	Yes
5	5.4	9.1.5	The studied indicator	ALT	30	Yes
	5.6		Dilution coefficient	1:50	30	No
6		9.1.6	Priority	S, R	1	Yes

Note.

In the sixth field, "Priority," "S" indicates severe cases and "R" (most commonly used) indicates standard cases.

The total number of fields in the protocol is 31. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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1.2.4.4 Recording the result (R)

Note.

The information in this section is provided in accordance with Section 10 of ASTM Protocol E1394. All test results should be presented using this type of the record.

Example:

R|1|^ALT^1:50|15.3|pg/mL|0 to 200|N|||||20180326172956<CR>

Table 1.2-8. Recording the result

Table 12: Comparison of Recording and Reporting						
Field		E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1		10.1.1	Message type	R	1	Yes
2		10.1.2	Serial number	1	5	Yes
3	3.4	10.1.3	The studied indicator	ALT	10	Yes
	3.6	10.1.3	Dilution coefficient	1:50	10	No
4		10.1.4	Result	115.3	12	No
5		10.1.5	Unit measured	pg/mL	10	No
6		10.1.6	Reference range	0 to 200	30	No
7		10.1.7	Result label	L,H,N,<, >	1	No
13		10.1.13	Time to complete the analysis	YYYYMMDD DHHMMSS	14	No

Note.

In the seventh field - "Result Label" - "L" applies if the result is below the normal range, "H" applies if the result is above the normal range, and "N" indicates that the result is within the normal range. "<" applies if the result is below the measurement range and ">" applies if the result is above the measurement range.

The total number of fields in the protocol is 14. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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1.2.4.5 Recording the query for information (Q)**Note.**

The information in this section is provided in accordance with Section 12 of ASTM Protocol E1394.

This record type is used to the query information about the indicators associated with a sample from the LIS system.

Example:

Q|1|^1234567||ALL|||||O<CR>

Table 1.2-9. Recording the query for information

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	12.1.1	Message type	Q	1	Yes
2	12.1.2	Serial number	1	6	Yes
3	12.1.3	Sample ID	^1234567	22	Yes
5	12.1.5	All indicators under study	All	10	Yes
13	10.1.13	Status of information query	O	1	Yes

Note.

The total number of fields in the protocol is 13. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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1.2.4.6 Finalizing the message (L)**Note.**

The information in this section is provided in accordance with Section 13 of ASTM Protocol E1394. Used in the last record of the transfer, indicates that the transfer has been completed.

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Table 1.2-11. The value of the symbols

Code	Description	Corresponding ASCII code
-->	Sending by software	
<--	Getting the software	

Where:

```
H|'^&||PSWD|MAGLUMI X3 |||||Lis||P|E1394-97|20100323<CR>
Q|1|^1234567||ALL|||||||O<CR>
L|1|IN<CR>
```

The software sends the above message to the LIS system using the ASTM E1394 data transfer protocol to request analysis information for the appropriate sample. The above message text performs the query for sample analysis information 1234567. This message is the opposite of the message in Section 1.2.5.2 of Appendix B.

which describes the process for receiving analysis information in response. Details on the meaning of the fields are provided in Appendix B, Section 1.2.4.

Note.

In the above message, the <ACK> character is an acknowledgement sent by the LIS system that requires the appropriate <ACK> command to be sent at the specified position. Otherwise, no connection to the LIS system will be registered in the software.

1.2.5.2 Response message with information about the analysis

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Table 1.2-11. The value of the symbols

Code	Description	Corresponding ASCII code
-->	Sending by software	
<--	Getting the software	

```

<--H|^&||PSWD| MAGLUMI X3 ||||Lis||P|E1394-97|20180319<CR>
<--P|1<CR>
<--O|1|1234567||^FT3|^FT4|R<CR>
<--L|1|N<CR>

```

Where:

The LIS system sends the above message using the ASTM E1394 data transfer protocol to send the analysis information of the corresponding sample in response. The text of the above message executes the sending of FT3 and FT4 analysis details for sample 1234567 in response. This message is the opposite of the message in Section 1.2.5.1 of Appendix B, which describes the process for requesting analysis details. Details on the meaning of the

fields are provided in Appendix B, Section 1.2.4.

Note.

In the above message, the <ACK> character is an acknowledgement sent by the analyzer that requires the appropriate <ACK> command to be sent at the specified position. Otherwise, no connection to the software will be registered in the LIS system.

1.2.5.3. Sending results

Content of the message

```

-><ENQ>
<--<ACK>
-><STX>
<--<ACK>
->H|^&||PSWD| MAGLUMI X3 ||||Lis||P|E1394-97|20180326<CR>
->P|1<CR>
->O|1|1234567||^FT3 ^FT4<CR>
->R|1|^FT3|123|pg/mL|0-200|N|||||20180326172956<CR>
->R|2|^FT4|25.1|pg/mL|0-50|N|||||20180326172956<CR>
->L|1|N<CR>
-><ETX>
<--<ACK>
-><EOT>
<--<ACK>

```

Table 1.2-11. The value of the symbols

Code	Description	Corresponding ASCII code
->	Sending by software	
<--	Getting the software	

Where:

```

H|^&||PSWD| MAGLUMI X3 ||||Lis||P|E1394-97|20180326<CR>
P|1<CR>
O|1|1234567||^FT3 ^FT4<CR>
R|1|^FT3|123|pg/mL|0-200|N|||||20180326172956<CR>
R|2|^FT4|25.1|pg/mL|0-50|N|||||20180326172956<CR>

```

The software sends the above message using the ASTM E1394 data transfer protocol to send the analysis results to the LIS system. The above text is the FT3 and FT4 results message for sample 1234567 that the analyzer sends to the LIS system. Details on the meaning of the fields are provided in Appendix B, Section 1.2.4.

1.3. HL7 Protocol

1.3.1. Description of the standard data transfer format

Format

<SB>CONTENT<EB><CR>.

Example

<SB>TEXT<EB><CR>

Table 1.3-1. Description of the standard data transfer format

Code	Description	Corresponding ASCII code
<SB>	The initial character	0x0B
T	The letter T	0x54
E	The letter E	0x45
S	The letter S	0x53
T	The letter T	0x54
<EB>	End symbol	0x1C
<CR>	Carriage return symbol	0x0D

1.3.2 Description of the separators

Table 1.3-2. Description of the separators

Separators	Description	Corresponding ASCII code
	Field separator	0x7C
—	Repeat separator	0x7E
^	Component separator	0x5E
&	The sub-component separator	0x26
\	Shielding symbol	0x5C

1.3.3 Description of the message types

Table 1.3-3. Description of the message types

Message label	Description
OML^O33	Sending analysis information in response
TSREQ	Query of analysis details
OUL^R22	Sending the analysis results
RSREQ	Query of analysis results

1.3.3.1 Response message with information about the analysis (OML^O33)

Note.

This message is used if the LIS system initiates sending orders to the software or if the LIS system sends a response to software orders.

Table 1.3-4. The structure of the OML^O33 message

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**Note.**

1. The total number of message segments in the protocol is 7. The list contains only the message segments used in this software. Other message segments can be added directly during data transmission; they will be automatically recognized by the software.

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1.3.3.3 Transfer Results (OUL^R22)

Note.

This message is used if the software initiates the sending of analysis results to the LIS system or if the software sends a response message with the results of the request to the LIS system.

Table 1.3-6. The structure of the OUL^R22 message

Message Segment	Value	Comment
MSH	The title of the message	
SPM	Sample	
ORC	General order segment	
{OBR	Research order	The recurring element
OBX	Study/result	The recurring element
[NTE]}	Comments	Optional, a repeating element

**Note.**

1. The total number of message segments in the protocol is 11. The list contains only the message segments used in this software. Other message segments can be added directly during data transmission; they will be automatically recognized by the software.

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1.3.4 Description of message segments

Note.

The total number of message segments in the protocol is 13. The list contains only the message segments used in this software. Other message segments can be added directly during data transmission; they will be automatically recognized by the software.

1.3.4.1 Segment Header (MSH)

Format:

MSH|^~\&|Sender||Recipient||The time of sending the message||Message type|Confirmation code|P|2.5|| NE|N E||UTF-8

Example:

MSH|^~\&|MAGLUMI||LIS||20151025105016||TSREQ|013329845|P|2.5||NE|NE||UTF-8

Table 1.3-8 Description of the MSH message segment

Position	Field	Content	Maximum length	Obliga tory	Comment
1	Field separator		1	Yes	Fixed value
2	Character encoding	^~\&	4	Yes	Fixed value
3	Sender application	MAGLUMI	180	Yes	
5	Recipient application	LIS	180	Yes	
7	The time of sending the message	20151025105016	14	Yes	YYYYMMDDHHMMSS
9	Message type	TSREQ	7	Yes	OML^O33, OUL^R22, RSREQ or TSREQ

10	Control element ID/Message confirmation code	013329845	20	Yes	Unique ID
11	Processing ID	P	3	Yes	Fixed value
12	ID version	2.5	8	Yes	Fixed value
15	Confirmation of receipt of the message	n/a	2	Yes	Fixed value
16	Confirmation from the application	n/a	2	Yes	Fixed value
18	Symbol set	UTF-8	16	Yes	Fixed value

**Note.**

The total number of fields in the protocol is 18. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.4.2 Patient Identification Segment (PID)

Format:

PID|1||Patient's ID||Patient's full name||yyyyMMdd(Date of birth)|Gender (Enumeration)|||||||||Age^The unit of measurement of age

Example:

PID|1||10248734||Zhang Jianguo||20110101|1|||||||8^1

Table 1.3-9 Description of the PID message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	No	1	1	Yes	Fixed value
3	Patient's ID	10248734	20		
5	Patient Name	Zhang Jianguo	50	Yes	
7	Date or time of birth	20110101	8		YYYYMMDD
8	Sex	1	1	Yes	1: male, 2: female, 3: unknown
38	38.1	Age	8	Yes	Positive integer
	38.2	Units of age measurement	1	Yes	1: year, 2: month 3: day, 4: hour

**Note.**

The total number of fields in the protocol is 38. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.4.3 Sample Details Segment (SPM)**Format:**

SPM|**#**|ID of sample ^ ^ID of control material^ID of sample rack^Sample rack position||**Sample type**|||||**Sample role**|||||

Example:

(sample) SPM|1|034581229010^101^A001^2||0|||||P|||||

(control material) SPM|1|^#2020030301#^C01^4||0|||||Q|||||

Table 1.3-10. Description of the SPM message segment

Position	Field	Content	Maximum length	Obligatoriness	Comment
1	Nº	1	1	Yes	Fixed value
2	2.1	Sample ID	034581229010	22	Yes
	2.3	ID of the control material	#2020030301#	22	Yes
	2.4	ID of the sample rack	A001/C01	22	
	2.5	Position of the sample rack	2/4	5	
4	Sample type	1/6	1	Yes	0: serum 1: Plasma 2: urine 3: hydrothorax and ascites 4: cerebrospinal fluid 5: Other
11	Role of the sample	P/Q	1	Yes	P: Patient sample Q: Control materials

**Note.**

The total number of fields in the protocol is 18. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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ORC|NW|||||0||201909131511433|||

Table 1.3-11. Description of the ORC message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	Status of the order	NW	2	Yes	NW: new physician's order RF: repeat order query
7	Priority	0	1	Yes	0: General 1: Quick check
9	Time of the order	201909131511433	14	Yes	yyMMddHHmmSS

Note.

The total number of fields in the protocol is 12. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

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Example: OBR|1||||FT3^1:10

Table 1.3-12. Description of the OBR message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	№	1	4	Yes	Starts at 1
4	4.1	Analysis code	FT3	Yes	Reagent channel no.
	4.2	Dilution coefficient	1:10		

Note.

The total number of fields in the protocol is 4. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

Configure the dilution factor sending function in the user software.

1.3.4.6 Result Information Segment (OBX)

Format:

OBX|||||Analysis code||Analysis result||Units of measurement||Reference range||Anomaly label||F||Begin time of analysis|| ||Module name||Equipment name||Analysis end time||Reference label||Message exchange type

Example:

OBX|1||FT4||20.9|mmol/L|8.9-17.2|H||F||20191025104016|||Module1^Maglumi|20191025104816|0|0

Table 1.3-13. Description of the OBX message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	№	1	4	Yes	Starts at 1
3	Analysis code	FT4	10	Yes	Abbreviated name of the reagent
5.1	Quantitative result	20.9	10	Yes	
6	Unit measured	2	10		Units of measure for reagent analysis results
7	Reference range	8,9–17,2	100		Reference range of values for reagent assay results
8	Abnormal value label	H	6		H: above normal HH: above the measuring range L: below normal LL: below the measuring range N: no label
11	Result Status	F	1	Yes	Fixed value, indicates that the result is final
14	Analysis start time	20191025104016	14	Yes	YYYYMMDDHHMMSS
18	18.1	Module Name	Module1	Yes	
	18.2	Equipment name	MAGLUMI	Yes	
19	Analysis end time	20191025104816	14	Yes	YYYYMMDDHHMMSS
20	Retest label	0	1	Yes	0: Result of standard analysis 1: Result of retest
21	Type of messaging	0	1	Yes	0: Automatic messaging 1: Manual messaging

**Note.**

The total number of fields in the protocol is 21. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.4.7 Segment with comments (NTE)

Format:

NTE|№

Example:

NTE|1

Table 1.3-14. Description of the NTE message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	№	1	4	Yes	Starts at 1

**Note.**

The total number of fields in the protocol is 3. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.4.8 Request Parameter Details (QPD) Segment

Format:

QPD|Query Name|||Sample ID^^Sample Rack ID^Sample Rack Position|Priority

Example:

QPD|RSREQ||123456^^N0001^2|0

Table 1.3-15. Description of the QPD message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	Name of the query message	RSREQ	7	Yes	RSREQ or TSREQ
3	3.1	Sample ID	123456	Yes	
	3.3	ID of the sample rack	N0001		
	3.4	Position of the sample rack	2		
4	Priority	0		Yes	0: General 1: STAT

**Note.**

The total number of fields in the protocol is 4. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.4.9 Response Control Details (RCP) Segment**Format:**

RCP|Query priority||| Response form

Example:

RCP|||R|

Table 1.3-16. Description of the RCP message segment

Position	Field	Content	Maximum length	Obligatory	Comment
1	Query priority	I	4	Yes	Fixed value, indicates the urgency of the query
3	Response Form	R	10	Yes	Fixed value, indicates that responses are given in real time

**Note.**

The total number of fields in the protocol is 3. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.

1.3.5 Description of the examples**1.3.5.1 Response message with information about the analysis****Actual message:**

```
<- MSH|^~&|Lis||MAGLUMI X3||20200311140350||OML^O33|2|P|2.5|||NE|NE||UTF-8
<- PID|1||34||zhangsan||19960101|Male|||||||||||||||||^Years Old
<- SPM|1||1234567^^^|||||P
<- ORC|NW|||||N
<- OBR|1|||FT3^
<- OBR|2|||FT4^
```

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Code	Description
<-	Messages received by this software

The software sends the above message using the HL7 data transfer protocol to send analysis information in response. The text of the above message executes the sending of FT3 and FT4 analysis details for sample 1234567 in response. This is an addendum to the message in Section 1.3.5.2 of Appendix B, which describes the process for requesting analysis details. Details on the meaning of the fields are provided in Appendix B, Section 1.3.4.

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Code	Description
-->	Messages sent by this software

The software sends the above message to the LIS system using the HL7 data transfer protocol to request the analysis details of the relevant sample. The above message text performs the query for sample analysis information 1234567. Details on the meaning of the fields are provided in Appendix B, Section 1.3.4.

Note.

After the software sends this message, the LIS system sends a reply message with the analysis information. See example in Appendix B, Section 1.3.5.1.

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Code	Description
-->	Messages sent by this software

The software sends the above message using the HL7 data transfer protocol to send the analysis results to the LIS system. The above text is the FT4 analysis results message for sample 1234567 that the analyzer sends to the LIS system. Details on the meaning of the fields are provided in Appendix B, Section 1.3.4.

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Code	Description
<--	Messages received by this software

The LIS system sends the above message to the software using the HL7 data transfer protocol to request the analysis result. The text of the above message is querying the result of sample 1234567. Details on the meaning of the fields are provided in Appendix B, Section 1.3.4.

**Note.**

After the software receives this message, it sends a reply message with the results of the analysis. See example in Appendix B, 1.3.5.3.

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Instruction for use

Shenzhen New Industries Biomedical Engineering Co., Ltd.

Dear user! Thank you for choosing the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories.

For safe and proper operation of the analyzer, as well as for effective use, please read the operating instructions carefully before starting to work with the device.

After reading this manual, keep it in an easily accessible place so that it can be accessed at any time.

If you have any questions about the Fully-auto Chemiluminescence Immunoassay Analyzer, contact your local representative office.

Shenzhen New Industries Biomedical Engineering Co., Ltd.

Address: No. 23, Jinxiu East Highway, Pingshan District, 518122 Shenzhen, People's Republic of China(RU)

Phone: +86-755-21536601 Fax: +86-755-28292740

Authorized Representative for the Russian Federation:

**LLC "KLS" (LLC CLS)125363, Moscow,
Novoposelkovaya str., 6 building 216.1 room 58
Phone: +7-495-9469597**

Shanghai International Holding Corp. GmbH (Europe)

Address: Elfestrasse 80, 20537 Hamburg, Germany

Phone: +49-40-2513175 Fax: +49-40-255726

Model	Catalog number
MAGLUMI X8	010101008801

Intellectual property rights

The intellectual property rights to this operating manual and the corresponding product belong to Shenzhen New Industries Biomedical Engineering Co., Ltd.

All rights reserved. Organizations and individuals are prohibited from copying, modifying or translating any part of this operating manual without the written permission of Shenzhen New Industries Biomedical Engineering Co., Ltd.

Snibe , **MAGLUMI Blossays** , **Biolumi** , **Preaccu** , **SnibeLinker** , **SnibeLis** , **SATLARS** , **Molecision**, etc . are registered trademarks or trademarks of Shenzhen New Industries Biomedical Engineering Co., Ltd.

DEVICE INFORMATION

Device Name: Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories

System device: the product consists of an chemiluminescence immunoassay module, a sample area, accessories and software. The chemiluminescence immunoassay module consists of a material supply module, a liquid circuit module, a temperature control module, a mechanical drive module, an optical signal registration module and an electrical circuit control module.

Appointment:

The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI® with accessories, in versions it is used with adapter reagents for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in blood serum, blood plasma, urine, whole blood.

GUARANTEE

The Analyzer manufacturer guarantees its compliance with the safety requirements for electronic and medical analyzers. The legislation establishes requirements for the installation and operation of the Analyzer, which must be followed by installation specialists and operators.

The manufacturer guarantees that the technical, functional and operational characteristics of the analyzers comply with the requirements of applicable standards and manufacturer's specifications during the life of the analyzer, provided that the consumer complies with the operating, transportation and storage conditions specified in the operational documentation.

Information about the User's Manual

Date of publication: 03/2023

Version: 1.1

Required software: ACCORDING to version 2.2.2.0 or other compatible version

Content

DEVICE INFORMATION	3
Warnings	1
INSTRUCTIONS FOR USE	2
DESCRIPTION OF WARNING SYMBOLS	2
PRECAUTIONARY MEASURES	2
PRECAUTIONS FOR USE	5
Chapter 1. Measurement principle	1
1.1 Stages of laboratory analysis	1
1.2 Measurement principle	2
1.3 Calibration method	3
Chapter 2. General information about the analyzer	1
2.1 General information.....	1
2.2 Specifications	1
2.3 General design of the analyzer	2
2.4 The main components of the analyzer	6
2.4.1 Chemiluminescence immunoassay module	6
2.4.2 SAMPLE AREA	21
2.5 Control unit.....	23
2.6 Labels	24
2.6.1 Module Labels	24
2.6.2 Other labels.....	30
2.7 The scope of delivery	32
2.8 Required materials (not included in the package)	32
Chapter 3. ANALYZER INSTALLATION	1
3.1 Storage and transportation requirements	1
3.1.1 Storage requirements	1
3.1.2 Transportation requirements.....	1
3.2 Installation Requirements	1
3.2.1 Requirements for spatial spacing.....	1
3.2.2 Environmental requirements.....	1
3.2.3 Power supply requirements	2
3.2.4 Water supply requirements	2
3.2.5 Drainage requirements.....	2
3.3 Unpacking and checking	3
3.3.1 Unpacking procedure.....	3
3.3.2 Movement and location	3
3.4 Analyzer Installation	3
3.4.1 Connection between modules	3
3.4.2 Electrical connection	4
3.4.3 Connecting washing liquid containers.....	4
3.4.4 Connection of liquid waste tanks.....	4
3.5 Starting and checking the system status	4
Chapter 4. BASIC OPERATING PROCEDURES	1

4.1	Standard operating procedures	1
4.2	Preparation for analysis	3
4.2.1	Pre-launch checks	3
4.2.2	Turning on and logging into the software.....	3
4.2.3	Replacement of consumables	3
4.2.4	Checking the analyzer status	10
4.2.5	Checking the conditions for performing the analysis	11
4.2.6	Loading of reagents	11
4.2.7	Uploading samples	14
4.3	Analyzing	15
4.3.1	Registration of calibration analysis	16
4.3.2	Registration of the QC analysis	17
4.3.3	Sample registration.....	17
4.3.4	Starting the analysis.....	18
4.3.5	Analysis of added samples	18
4.4	Tracking the execution of the analysis	19
4.5	Working with results	19
4.6	Completion of the analysis	20
4.6.1	Preparing to put the analyzer into sleep mode or turn it off	20
4.6.2	Power off.....	20
4.6.3	Actions after power failure.....	20
4.6.4	Removal of the reagent kit.....	21
4.6.5	Removing and maintaining the sample rack.....	21
4.6.6	Waste disposal.....	22
Chapter 5.	Software functions.....	1
5.1	Description of the software interface.....	1
5.1.1	Main menu.....	1
5.1.2	Interfaces	2
5.1.3	Submenu.....	2
5.1.4	Toolbar	3
5.1.5	Dialog boxes.....	5
5.1.6	Buttons.....	5
5.2	Function description [Sample].....	5
5.2.1	Sample registration	6
5.2.2	Scanning	8
5.2.3	Application for analyses	8
5.2.4	Worklist	9
5.3	Description of the screen functions [Result]	10
5.3.1	Search for results	11
5.3.2	Deleting the results	11
5.3.3	Information about the result	11
5.3.4	Repeat measurement.....	12
5.3.5	Sending results.....	12
5.3.6	Recalculation of results	12

5.3.7	Exporting results.....	13
5.3.8	Printing the results.....	13
5.3.9	Description of result labels.....	13
5.4	Function description of the [Reagent] screen.....	14
5.4.1	Information about the reagent.....	14
5.4.2	Curve query	17
5.4.3	Calibration registration.....	17
5.4.4	Calibration results.....	17
5.4.5	Information about the reagent.....	18
5.4.6	Reagent Information List.....	20
5.4.7	Reagent consumption statistics.....	21
5.4.8	Enabling and disabling the use of reagents.....	21
5.5	Description of the screen functions [QC]	22
5.5.1	QC registration	22
5.5.2	QC results	24
5.5.3	Analysis of control materials.....	24
5.6	Description of the screen functions [Utility]	25
5.6.1	System checks	25
5.6.2	Loading liquid to check the light signal	27
5.7	Description of the functions of the [Settings] screen.....	27
5.7.1	Setting Parameters	28
5.7.2	System Settings	39
5.7.3	System maintenance	44
5.8	Description of [Home screen] screen functions.....	48
5.8.1	Manual of operations.....	48
5.8.2	Module monitoring.....	52
5.9	Description of the functions of the [Status] screen.....	53
5.9.1	Status of the chemiluminescence immunoassay module.....	53
5.9.2	Status of the sample processing system.....	62
Chapter 6.	System maintenance	1
6.1	Tools and materials.....	1
6.1.1	Accessories (supplied with the analyzer)	1
6.1.2	Materials to be prepared	2
6.1.3	Cuvettes	2
6.1.4	Disposable tip	2
6.1.5	Washing liquid.....	2
6.1.6	Starting reagent.....	2
6.1.7	Solution for cleaning system tubes	3
6.1.8	Purified water	3
6.2	Daily maintenance.....	3
6.3	Weekly maintenance	3
6.4	Monthly maintenance	4
6.5	Prolonged non-use.....	4
Chapter 7.	The main malfunctions and ways to eliminate them	1

7.1	Classification of faults	1
7.2	Faults not accompanied by an alarm and how to rectify them	1
7.3	Faults not accompanied by an alarm and how to rectify them	13
Chapter 8.	Environmental requirements.....	13
8.1	Disposal and destruction.....	13
Chapter 9.	List of international standards applicable to MD.....	13
Chapter 10.	Storage, transportation, packaging	15
Appendix A.	Operational characteristics	1
1.1.	Accuracy and repeatability of dosing	1
1.2	Accuracy and deviations of the reaction disk temperature	1
1.3	LC check components	1
1.4	Transfer	2
1.5	Reproducibility of clinical assays within a single batch.....	2
Appendix B.	Network interfaces	1
B.1	Protocol	1
B.2	Protocol ASTM	1
B.2.1	Description of control codes.....	1
B.2.2	Description of the standard data transfer format	1
B.2.3	Description of the separators	2
B.2.4	Description of the message types	2
B.2.5	Example.....	8
B.3	HL7 Protocol	13
B.3.1	Description of the standard data transfer format	13
B.3.2	Description of the separators	13
B.3.3	Description of the message types	13
B.3.4	Description of message fields.....	15
B.3.5	Example	23

Warnings

This section contains important information on safety and proper operation of the device, as well as information on regulatory requirements.

Please read the operating instructions before using the analyzer.

INSTRUCTIONS FOR USE

The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI® with accessories, in versions it is used with adapter reagents for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in blood serum, blood plasma, urine, whole blood.

This user manual is intended for users of the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories. The manual describes the principles of the device, operation and operation of the analyzer, as well as provides information on routine maintenance procedures, troubleshooting measures and other aspects of using the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories. The analyzer must be operated in accordance with this manual.

Indications and contraindications for the use of the analyzer:

Indications:

Indications for the use of the analyzer MAGLUMI® are the need for comprehensive diagnosis and monitoring of specific diseases* and medical conditions.

It is used in clinical practice for qualitative and/or semi-quantitative and/or quantitative analysis of analytes in blood serum, blood plasma, urine, whole blood.

*Testing: including hormones, tumor antigens, proteins and peptides, liver diseases, cardiomyopathy, myocardial diseases, immune function, autoantibodies, infectious diseases, concentration of vitamins and drugs in the blood, other enzymes, D-dimer, coagulation, microalbuminuria, human immunoglobulin E, prolactin.

Contraindications:

Contraindications to the use of a medical device for in vitro diagnosis and monitoring have not been identified.

DESCRIPTION OF WARNING SYMBOLS

Symbol	Text	Description
	Warning	Warnings indicate the presence of a condition in which there is a risk of injury, damage to the analyzer, as well as data errors or biological hazards.
	Note	The note contains important information that you should pay attention to.

PRECAUTIONARY MEASURES

For safe operation of the analyzer, it is necessary to read the instructions carefully before starting work. Failure to follow these precautions may result in injury or damage to the analyzer.

This system is manufactured in accordance with the safety requirements for the production of electronic and medical analyzers. There are statutory requirements and standards for the installation and operation of this system, which must be followed by installation specialists and operators.

Warning

1. The user must regularly inspect and maintain the device in strict accordance with the instructions in order to avoid disruption of the analyzer and the occurrence of a threat to the safety of personnel.
2. To ensure the safe and reliable operation of the analyzer, its installation and maintenance must be carried out only by authorized technical personnel or with the permission of such specialists. All analyzer components must be purchased and tested by our company or authorized dealers.

1. Warning of injuries caused by moving parts

To avoid injury caused by moving parts during operation of the analyzer, observe the following precautions.

Warning

1. Do not touch the moving parts of the analyzer and do not interfere with their movement during the operation of the system. Do not place foreign objects in the way of moving parts. Failure to comply with this requirement may result in injury or damage to the analyzer.
2. To avoid injury or damage to the analyzer, when loading liquid and solid consumables into the analyzer, follow the established rules and properly remove the sealing film, lids or other packaging from the consumables.

2. Warning of the danger of electric shock

To avoid electric shock, observe the following precautions.

Warning

1. Make sure that the input voltage meets the requirements for the operation of the analyzer.
2. Only an authorized service technician may open the front or back cover and other elements during operation of the analyzer.
3. Leakage of reagents, samples and other liquids inside the analyzer may cause malfunction or electric shock. If liquid gets into the analyzer, immediately turn off the power supply of the device and contact the SNIBE maintenance specialists.
4. If you need to open the front or back cover of the analyzer to replace parts, turn off the analyzer power before opening.
5. Improper grounding may result in electric shock and damage to the analyzer.

3. Fire hazard warning

The use of organic solutions can lead to a fire. Observe the following precautions.

Warning

1. Do not use organic solutions for analysis.
2. The analyzer design is not explosion-proof. Do not use organic solutions with an ignition temperature of less than 65 °C near the analyzer. Otherwise, it may

cause a fire or explosion.

4. Warning of biochemical hazards

To avoid the occurrence of biochemical hazards, observe the following precautions.

Warning

1. Before using, carefully read the safety rules and precautions given in the instructions for reagents and consumables.
2. Improper handling of samples can lead to infection. It is forbidden to work with samples, mixtures and liquid waste with bare hands. Contact of these materials with other parts of the body is also not allowed. To avoid infection, always wear gloves, a gauze mask and special clothing while working; if necessary, use protective glasses.
3. The use of certain types of reagents and washing liquid can lead to skin lesions. When working with reagents and washing liquid, care should be taken to avoid their contact with the skin of hands and clothes. In case of accidental contact of reagents or washing liquid on clothes or skin of hands, wash them immediately with soap and water. In case of accidental contact with the eyes, rinse them immediately with plenty of water and consult an ophthalmologist.
4. If a small amount of reagent or sample gets on the analyzer during operation, use a clean cotton cloth soaked in alcohol to clean the device before starting work. If a large amount of reagent or sample gets on the analyzer, immediately stop using it and contact an engineer authorized by our company to promptly solve the problem.
5. Before transporting the analyzer over a long distance, it must be thoroughly disinfected to avoid the risk of infection.

5. Liquid waste disposal

To avoid environmental pollution as well as personal injury when disposing of liquid waste, observe the following precautions.

Warning

1. The disposal of certain substances contained in quality control and calibration solutions and liquid wastes is subject to contaminant regulations and disposal standards. Observe local waste disposal regulations and requirements specified by the manufacturer or distributor of the respective reagent.
2. When analyzing samples from patients with infectious diseases, liquid waste should be disposed of by the healthcare facility's infected waste treatment department.

6. Disposal of used analyzers

Observe the following requirements when disposing of used analyzers.

Warning

Disposal of certain substances contained in spent analyzers is subject to regulations regarding contaminants. Observe the relevant local regulations when disposing of

the product.

7. Protection from computer viruses

To avoid infecting your computer with viruses, observe the following precautions.

Warning

1. To prevent virus infection or damage to the computer software due to improper use or other reasons, use the data transfer function in accordance with the restrictions set by our company. Computer viruses can be transmitted via floppy disks, U-drives, networks, and other channels.
 2. Do not install any software or hardware other than those specified by our company, as this may result in incorrect operation of the software. Do not run other programs while the system is running.
-

PRECAUTIONS FOR USE

For proper and effective use of the analyzer, the following precautions should be read carefully.

1. General precautions

Before using this analyzer, its purpose and general precautions should be established. It is not possible to guarantee safe operation of the analyzer if the rules specified in this operating manual are not followed.

Note.

1. Clinical symptoms and test results from other laboratories should be considered when making clinical decisions based on test results.
 2. The User Manual is subject to change without prior notice. If necessary, you should contact a customer service representative.
 3. This device should only be used by trained personnel who have knowledge of basic analyzer operation, analysis and laboratory skills. The device is intended for use only under appropriate operating conditions.
 4. Do not touch your computer monitor, mouse or keyboard if your hands are contaminated with chemicals.
 5. Do not bend or squeeze the drain tubes. Failure to do so may cause them to become clogged and consequently drain liquid waste from other components, resulting in damage to the analyzer in case of serious irregularities.
 6. The analyzer generates heat during operation, which is expelled through holes in the rear panel of the analyzer. The workroom must be well ventilated to ensure heat dissipation. If necessary, a ventilating device should be used. However, the analyzer should not be exposed to direct air current. Failure to do so may result in inaccurate analytical results.
 7. Before starting the analyzer for the first time, adjust the settings of all components to ensure the accuracy of their parameters.
 8. Samples, reagents, starting reagents and flushing fluid must be free of air bubbles. Otherwise, this may lead to errors in sample addition and the validity
-

- of the analytical results is not guaranteed. Do not move or change reagents until reagent addition is complete.
9. Start reagent containers in the chemiluminescence immunoassay module must be properly connected. Do not use mixed starting reagents. Otherwise, the validity of the analytical results is not guaranteed.
 10. To ensure safe analyzer operation and accurate test results, use compatible cuvettes and disposable tips approved by our company.
 11. To ensure safe analyzer operation and consistent test results, do not use consumables that are past their expiration date.
 12. The analyzer should be turned on at least 30 minutes prior to operation to ensure stabilization of the measurement system.
 13. Ensure that all consumables (specifically reagents, washing fluid, cuvettes, and disposable tips) are adequate before performing any analysis. Check that the water supply to the analyzer meets the water supply requirements specified in Section 3.2.4.
 14. Before performing any analysis, perform a background wash (BGW) at least once and ensure that the results obtained are within the normal range. Otherwise, the validity of the analytical results is not guaranteed.
 15. Quality control procedures must be followed prior to analyzing the sample. Otherwise, the validity of the analytical results is not guaranteed.
 16. The analyzer is equipped with an alarm indicator light and speaker. If an analyzer error occurs, the alarm light and speaker will sound an appropriate alarm. To ensure proper analyzer operation and accurate test results, take appropriate corrective action immediately after an error occurs. For troubleshooting information, refer to Chapter 7 of this manual: "Basic Troubleshooting and ways to fix them".
 17. Use only SNIBE-approved tubing system flushing fluid to clean and maintain the analyzer tubing.
 18. In the event of spillage or ingress of hazardous substances on the surface or in the analyzer housing, properly disinfect the surfaces involved. Do not use cleaners or disinfectants that may react chemically with component materials or parts of the analyzer. If you have any questions regarding the compatibility of a disinfectant or cleaner with the analyzer materials or components, contact the manufacturer or their representative.
 19. The technicians authorized by our company carry out the calibration in accordance with the established rules.
 20. In the event of a serious incident involving the use of the device, notify Shenzhen New Industries Biomedical Engineering Co., Ltd. (Snibe) or our authorized representative and the competent authority of the Member State in which you are located.
 21. If the analyzer is broken or needs to be disposed of, the software of the Fully-auto chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories must be uninstalled and the data purged.
-

2. Operational conditions

Note.

When installing the analyzer, observe the workspace requirements specified in the User's Manual. Installation and use of the analyzer under conditions other than those specified may result in unreliable results or damage to the analyzer.

3. System maintenance

Warning

1. Maintenance of the analyzer should be performed on a regular basis following the requirements specified in the appropriate sections of this manual. Improper maintenance can adversely affect the accuracy and reproducibility of analytical results and may result in analyzer malfunction or injury.
2. The surface of the analyzer may become dusty if the analyzer is not used for a long period of time. To clean the surface, dampen a clean cotton cloth with water, wring it out, and gently wipe the surface of the analyzer. Take the necessary precautions to prevent water from entering the analyzer housing.
3. Disconnect power to the analyzer before performing maintenance (unless performing live maintenance). Failure to do so may result in analyzer failure or injury.
4. The analyzer may receive samples from potentially infected patients. Always wear gloves and special clothing during maintenance to prevent infection.
5. The analyzer does not contain any user-serviceable components. Do not attempt to open the case of the device or remove any components. If necessary, contact an authorized SNIBE service technician for assistance.
6. Contact an authorized service technician to check if the overheat protection system, the fluid level sensor protector, or the entire system is working.

4. Samples, reagents and quality control solution

Warning

1. The presence of substances such as drugs, anticoagulants, and preservatives in the sample may affect the results of the analysis.
2. The presence of hemolysis, lipemia, and jaundice affects test results.
3. Make sure there are no clots in the sample. Otherwise, the disposable tips may become clogged, significantly affecting the results of the assay.
4. Samples should be stored properly. Improper storage of samples can lead to changes in their component composition and incorrect analytical results.
5. To prevent samples from evaporating, do not leave them open for long periods of time. Evaporation of the sample may result in incorrect analytical results.
6. Improper storage of reagents and quality control solutions can result in incorrect test results and improper analyzer operation, even if these solutions are not expired. Quality control reagents and solutions should be stored and used in accordance with the instructions of the respective reagent

manufacturer.

7. After changing the reagent, calibration must be performed. Lack of calibration can adversely affect the accuracy of analytical results.
-

5. Data backup

Note.

The system has a function to automatically save data to the computer's hard disk. However, if data is deleted or the hard disk is damaged, the data stored on the disk cannot be recovered. Therefore, you should regularly back up the analytical results and analyzer parameters and save them on other media such as a CD-ROM drive.

6. Electromagnetic compatibility

Note.

1. The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories meets the requirements for electromagnetic radiation levels and immunity to interference specified in IEC 61326-2-6 and IEC 61326-1.
 2. For proper operation of the analyzer, it is the responsibility of the user to ensure that the analyzer is operated in accordance with EMC requirements.
 3. It is recommended to evaluate the electromagnetic environment before using the analyzer.
-

Warning


1. The analyzer is designed and tested in accordance with the CISPR 11 standard for Class A devices. It is not intended for residential use and may not provide an adequate level of protection against radio emissions in such environments.
 2. This equipment is intended for use in a professional medical facility. In a home use environment, the device will most likely not work properly. If electromagnetic interference is suspected, proper operation can be restored by increasing the distance between the equipment and the source of the interference.
 3. Do not use this unit close to sources of strong electromagnetic radiation (such as sources of intentional RF radiation), as this may affect the correct operation of the device.
-

CHAPTER 1 Measurement principle


The operating principles of the analyzer include the steps of laboratory analysis, the principle of measurement, and the calibration method. This section provides a detailed description of these procedures.

1.1 Stages of laboratory analysis

One step assay

Notes

Two step assay


Fig. 1.1-1. Stages of laboratory analysis

1.2 Measurement principle

The analyzer, whose operation is based on the principles of immunoreaction, uses heterogeneous magnetic separation technology combined with chemiluminescent reaction to directly label antigens or antibodies with ABEI luminescent tags. This allows to induce a reaction with antigens or antibodies in the analyzed sample so that these antigens or antibodies immobilized by the magnetic particles form an antigen-antibody immune complex. After the free antigen or antibody and other components are removed in the magnetic separation process, a starting reagent is added to the antigen-antibody immune complex to make the tag fluorescent. The luminescence level is then measured from the ROC curve to determine the concentration of the analyte.

The analyzer's photomultiplier tube (PMT) is used to record the light flux emitted during the chemiluminescence reaction. The linear measurement range of the light wavelength recorded by the FEU is 300-650 nm, and the maximum light wavelength emitted from the chemiluminescent reaction is 420 nm. The light signal generated during the chemiluminescence reaction enters the PMT, where it is converted into an electrical signal. Peripheral circuits collect, process, and convert this signal into a digital value representing the magnitude of the light intensity. The average life of the photomultiplier tube of this product is 27,000 hours.

To eliminate discrepancies between different PMTs and to ensure uniformity of results obtained with different analyzers, relative light units (RLU) are used as the unit of measurement of the raw data.

The following describes the process for performing the analysis in the analyzer. After the sample and reagent are added to the cuvette, they undergo three stages of incubation, washing, and separation before the cuvette is placed in the chamber. Starting reagent 1 is injected into the cuvette, after 2.5 s starting reagent 2 is injected into the same cuvette to start the chemiluminescence reaction. The recording of optical signals by the device begins 0.1 second after the chemiluminescent reaction and lasts for 3.0 seconds, as shown in Figure 1.2-1. The measurement time depends on the incubation time and the order in which the reagent is handled. If the one-step method is used and the incubation time is 1 minute, the measurement time (the period required to obtain the first result) is no more than 7 minutes. If the two-step method is used and the incubation time is 5 minutes, the measurement time (the period required to obtain the first result) will be no more than 11 minutes.


Fig. 1.2-1. Kinetic curve of chemiluminescence reaction

1.3 Calibration method

Since the actual operating conditions of the analyzer differ from laboratory conditions, it is necessary to adjust the reference curve in order to obtain a working curve corresponding to the actual operating conditions.

Short description:

- The reference curve is calculated based on 10 standard points.
- The values of the relative light units at the two calibration points obtained as a result of calibration are compared with the values of the relative light units of the corresponding concentration on the reference curve.
- The difference between the relative light intensity at the two calibration points and the relative light intensity of the corresponding concentration on the reference curve is then calculated. After that, a linear relationship is established based on the recalculated relative light intensity (Y axis) and concentration (X axis).
- For other points of the reference curve, the difference of the corresponding relative light units is calculated using a compensation curve, after which the values of relative light intensity (Y axis) and concentration (X axis) are recalculated.
- The recalculated curve is used as a valid working curve (see Fig. 1.3-1. "The principle of calibration").


Fig. 1.3-1. The principle of calibration

CHAPTER 2 General information about the analyzer

2.1 General information

The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories and a set of compatible diagnostic reagents represent an accurate system that allows direct immunochemiluminescent analysis based on the separation of compounds labeled with aminobutyl-ethyl-isoluminol (ABEI) in a magnetic field. It is used to carry out qualitative and (or) quantitative analysis of the content of analytes in samples of human origin. The analyzer automatically performs sample dosing, reagent loading, incubation, rinsing, measurements and calculation of results, which reduces the number of errors in the analysis results, as well as increases their accuracy and reproducibility.

2.2 Specifications

Table. 2.2-1. Analyzer Specifications

Component		Specifications
Main characteristics	Test execution speed	600 tests per hour
	Sample type	Serum, plasma, urine, whole blood
	Barcode Type	Code 128, Code 39, Code 93, Codabar, 2/5 Interleaved
	QC	Monthly quality control
Device Characteristics	Sample volume	10-280 μ l
	Sample area	30 positions for sample stands, a total of 300 positions for sample tubes
	The sampler	Liquid level detection, liquid level tracking, detection of clots, bubbles and collisions
	Reagent Dispenser Holder	Automatic washing, liquid level detection, liquid level tracking, collision detection
	Reagent area	There are 42 positions for reagents in total
	Reagent volume	10-450 μ l
	Temperature in the reagent area	Operating temperature: 8-12 $^{\circ}$ C, storage temperature: 2-8 $^{\circ}$ C
	Number of cuvettes	The total number of loaded cuvettes is 2,912
	Number of disposable tips	The total number of loaded cuvettes is 576

General information about the analyzer

	Reaction temperature	37.0 °C (± 0.3 °C), the deviation does not exceed 0.1 °C
	Mixing method	Mechanical vibrations
External interfaces	Interface	Ethernet
Operating conditions	Temperature	10–30 °C
	Relative humidity	≤ 70 %
	Other	Do not place the device near sources of electromagnetic interference
Storage conditions	Temperature	-20 to 55 °C
	Relative humidity	≤ 93 %
	Atmospheric pressure	50.0–106.0 kPa
	Other	Bright sunlight and corrosive gases should not enter the room, It is necessary to ensure good ventilation
Safety category	Degree of protection against electric shock	Type I
	Type of overvoltage	Type II
	Degree of contamination	Degree 2
Analyzer	Power	2200 V·A
	Power supply	100-240 V AC, 50/60 Hz
	Weight	670kg
	External dimensions (L*W*H)	192 x 118 x 150 cm
	Package Dimensions (L*W*H)	Chemiluminescence immunoassay module: 190 x 140 x 185 cm Sample area: 140 x 90 x 129.5cm
	Unless otherwise specified, the size and weight tolerances are $\pm 5\%$.	

2.3 General design of the analyzer

The device consists of an chemiluminescence immunoassay module, a sample area, accessories and software. The chemiluminescence immunoassay module consists of a material supply module, a liquid module, a temperature control module, a mechanical drive module, an optical signal registration module and an electrical circuit control module.

- The material supply module consists of a compartment for storing cuvettes, a compartment for storing tips, a reagent area, a tray for washing liquid, a container for starting reagents, a liquid

waste module, a container for used cuvettes and a container for used tips.

- The liquid module consists of the liquid contour of the sample dispenser, the liquid contour of the reagent dispenser, the liquid contour of the flushing unit, the liquid contour of the chamber and the liquid contour of the reagent area.
- The temperature control module consists of a heating module for the reaction disk, a cooling module for the reagent area and a temperature control module for the photoelectronic multiplier.
- The mechanical drive module consists of a device for transporting the sample area, a sample conveyor of the SPS system, a compartment for storing cuvettes, a compartment for storing tips, a reaction disk, a sample disk, a sampler, a reagent probe, a flushing unit, a chamber, a capture device for new cuvettes, a capture device for samples, an incubator capture device, a flushing capture device a node, a chamber capture device, and a shaker of the reagent area.
- The optical signal registration module consists of a chamber module, a PMT module, and a main control loop module.
- The electrical circuit control module consists of a power supply module, a main control board, a cable system, a number of sensors, and motors.


Fig. 2.3-1. Analyzer MAGLUMI X8: side view

(1) The trial platform

(2)
Chemiluminescence
immunoassay
module


Fig. 2.3-2. Analyzer MAGLUMI X8: Front view


Fig. 2.3-3. Analyzer MAGLUMI X8: Rear view

(1) High concentration waste

(2) Pure water input of the

(3) A hole for the removal

outlet of the chemiluminescence
immunoassay module

chemiluminescence immunoassay
module

of low-concentration waste
from the
chemiluminescence
immunoassay module

(4) Power switch of the main
installation of the
chemiluminescence immunoassay
module

(5) Power switch of the
chemiluminescence
immunoassay module
conveyor

2.4 The main components of the analyzer

A detailed list of the individual components of each module of the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories is presented below.

2.4.1 Chemiluminescence immunoassay module


Fig. 2.4-1. Chemiluminescence immunoassay module: side view


Fig. 2.4-2. Chemiluminescence immunoassay module: front view

- | | | |
|---|--------------------------|--------------------------|
| 1) Indicator of the quantity of consumables | 2) Alarm indicator light | 3) Emergency stop button |
| 4) Starting reagent 2 | 5) Starting reagent 1 | 6) Washing liquid |
| 7) RFID tag reader | | |


Fig. 2.4-3. Chemiluminescence immunoassay module: rear view

(1) High concentration liquid waste outlet

(2) Thermal air outlet of the cooling module of the reagent area

(3) Reserved port

(4) Low concentration liquid waste outlet

(5) Power switch, power interface and network interface

(6) Thermal air outlet of the power supply


Fig. 2.4-4. The internal structure of the chemiluminescence immunoassay module

- | | | |
|--|---|---|
| (1) Tips storage compartment | (2) Cuvette storage compartment | (3) SPS system Sample Conveyor |
| (4) Sample disk | (5) Sampler holder | (6) Sample capture device |
| (7) Chamber | (8) Chamber capture device | (9) Flushing unit |
| (10) Pump | (11) The device for capturing the block of tips | (12) New cup capture device |
| (13) Incubator capture device | (14) Reaction disk | (15) Holder of the left reagent dispenser |
| (16) The left area of the reagents | (17) Radio Frequency Tag Reader | (18) The device for capturing the flushing unit |
| (19) Holder of the right reagent dispenser | (20) The right area of the reagents | (21) Footing |
| (22) Sample area | | |

2.4.1.1 Reaction disk

Warning

Do not touch the reaction disk while the analyzer is running. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-5. Schematic drawing of the reaction disk

- | | | | |
|--|--|--|---------------------------------|
| 1. The inner ring of the reaction disk | 2. The outer ring of the reaction disk | 3. Insulation cover of the reaction disk | 4. Stirrer of the reaction disk |
|--|--|--|---------------------------------|

Function

The outer ring of the reaction disk: moving the cuvette to different working areas — to the position of introducing the reagent and mixing, to the position for transportation to the inner and outer ring, etc.

Inner ring of the reaction disk: incubation of a mixture of sample and reagent in a cuvette.

Reaction disc agitator: mixing of components in cuvettes into which samples and reagents have already been added.

2.4.1.2 Reagent area

Components


Fig. 2.4-6. Schematic drawing of the reagent area

(1) Radio Frequency
Tag Reader

(2) Cover

(3) Place for reagents

(4) Shaker

Function

The reagent area is designed to accommodate a set of reagents; the reagent dispenser aspirates reagents through holes on the lid. When the analyzer is in standby mode, the reagent area cover is closed. During the operation of the analyzer, the lid of the reagent area opens automatically. If it is necessary to turn off the cooling system, first turn off the chemiluminescence immunoassay module.

2.4.1.3 Reagent Dispenser Holder

Warning

It is forbidden to touch the holder of the left and right reagent dispensers during the operation of the analyzer. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-7. Schematic drawing of the left reagent dispenser

(1) X Axis

(2) Y axis

(3) Z axis

(4) Reagent probe

The device and principles of operation of the holder of the left and right reagent dispenser are similar. As an example, only the holder of the left reagent dispenser is considered in detail in this manual.

Function

The reagent probe moves reagents from the reagent kit into cuvettes in the reagent area. The following functions are available: automatic washing, liquid level detection, liquid level tracking and collision detection.

2.4.1.4 Sampler

**Warning**

It is forbidden to touch the mechanism of the sampler holder during the operation of the analyzer. Otherwise, injury or damage to the analyzer may occur.


Components

Fig. 2.4-8. Schematic drawing of the sampler

(1) Liquid displacement module

(2) Rotating holder

(3) X Axis

Function

The sampler performs loading of disposable tips, aspiration of a fixed amount of sample, calibration solution, quality control solution, etc. from the sample conveyor of the SPS system or the sample disk, and then dispenses them into the cuvettes. It is equipped with the functions of liquid level detection, liquid level tracking, detection of clots, bubbles and collisions, as well as other functions.

2.4.1.5 Sample disk

Warning

Do not touch the sample disc while the analyzer is running. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-9. Schematic drawing of the sample disk

(1) Support

(2) Drive mechanism

(3) Rotating sample disc platform

Function

The sample disc moves the cuvettes to various areas where the sampler holder dispenses samples, or to other work areas, and also serves as a buffer zone for cuvettes with already added samples.

2.4.1.6 Storage compartment for the cuvette

Warning

Do not touch the storage compartments and the upper return compartment in the cuvette storage compartment while the analyzer is in operation. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-10. Schematic drawing of the cuvette storage compartment

- | | | |
|---|---------------------|--------------------------|
| (1) Sliding compartment | (2) Support panel | (3) Reagent storage area |
| (4) The area of use of the new cuvette unit | (5) The return area | |

Function

The cuvette storage compartment allows you to move new cuvette blocks to the area for use, as well as return used cuvettes.

This compartment consists of three parts: the lower and upper compartments of the sliding compartment and the area for using new cuvette blocks located above the upper compartment for returning blocks. Nine cuvette blocks can be stored in the lower compartment, and six cuvette blocks in the upper compartment. One cell block can be placed in the area of use of the new cell block. Thus, in general, it is possible to store 16 blocks of cuvettes. Each block contains 182 cuvettes. A total of 2,912 cuvettes can be loaded into the chemiluminescence immunoassay module.

The compartment for the return of the cuvette consists of two parts: the lower and upper compartments. Nine spent cuvette blocks can be stored in the lower compartment of the return

compartment, and seven spent cuvette blocks can be stored in the upper compartment. Thus, a total of 16 spent cell blocks can be stored in the chemiluminescence immunoassay module.

Five blocks of new, unused or spent cuvettes can be placed in the buffer container of the retractable compartment of the chemiluminescence immunoassay module.

2.4.1.7 Tips storage compartment

Warning

It is forbidden to touch the tip block gripper and the transfer mechanism while the analyzer is operating. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-11. Schematic drawing of the tip storage compartment

(1) Sliding compartment

(2) The device for capturing the block of tips

(3) Buffer container

(4) Transfer mechanism

Function

The tip storage compartment allows you to move blocks with disposable tips to the area for use, as well as return used tip blocks. The retractable tip storage compartment can hold four blocks of disposable tips. One block of disposable tips can be stored in the buffer container and the transfer mechanism. There

are 96 disposable tips in one unit. A total of 576 disposable tips can be loaded into the chemiluminescence immunoassay module.

2.4.1.8 Flushing unit

Warning

Do not touch the flushing unit while the analyzer is in operation. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-12. Schematic drawing of the flushing unit

- | | | | |
|-------------------------------|---------------------------|-----------------------------|----------------------------------|
| (1) Disc of the flushing unit | (2) Primary flushing unit | (3) Secondary flushing unit | (4) Agitator of the washing unit |
|-------------------------------|---------------------------|-----------------------------|----------------------------------|

Function

The washing unit is primarily designed to remove components of the reaction solution that are not associated with a labeled antigen or antibody of magnetic microparticles in the solid phase by means of three-time washing. The disk of the flushing unit moves the cuvettes into the magnetic field and other working areas. The first and second needles for removing liquid waste from the primary flushing unit, as well as the third and fourth needles for removing liquid waste from the secondary flushing unit are designed to aspirate flushing liquids and drain them in the flushing unit. Four flushing needles are designed to inject the flushing liquid. The agitator of the washing unit is designed to mix the mixture in a cuvette.

2.4.1.9 Chamber

Warning

Do not touch the chamber components while the analyzer is running. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-13. Schematic drawing of the chamber

(1) Liquid waste disposal mechanism

(2) Tip for starting reagent 1

(3) Tip for starting reagent 2

(4) The PMT module

(5) Chamber disc

Function

In the chamber, the starting reagent 1 and 2 are introduced into a cuvette that has passed the washing cycle, then an optical signal is generated as a result of their action on the analyte. This signal is detected, processed and converted into a value representing the light intensity by means of a PMT and a control loop. The liquid waste disposal mechanism is designed to remove liquid waste from cuvettes for which the analysis procedure has been completed.

2.4.1.10 Cuvette capture device

Warning

It is forbidden to touch the cuvette capture devices while the analyzer is running. Otherwise, injury or damage to the analyzer may occur.

Components

The cuvette capture device consists of five independent components: a capture device for new cuvettes, a capture device for samples, an incubator capture device, a flushing unit capture device and a chamber capture device. The design and principles of operation of these components are similar. As an example, this manual discusses in detail only the gripper device for new cuvettes.


Fig. 2.4-14. Schematic drawing of the gripper device for new cuvettes

(1) X Axis

(2) Y axis

(3) Z axis

(4) Clamping mechanism

Function

The gripper device for new cuvettes is primarily designed to move the cuvettes from the cuvette storage compartment to the sample disk.

The sample capture device is primarily designed to move the cuvettes from the sample disk to the outer ring of the reaction disk.

The incubator capture device is primarily designed to move the cuvettes from the outer ring of the reaction disk to the inner ring of the reaction disk.

The capture device of the flushing unit is primarily designed to move the cuvettes from the inner ring of the reaction disk to the flushing unit or move the cuvettes from the flushing unit to the outer ring of the reaction disk.

The chamber capture device is primarily designed to move the cuvettes from the flushing unit into the

chamber.

The movement of the cuvettes between the various components of the analyzer is achieved through the interaction of various capture devices. All capture devices have the function of dumping spent cuvettes into the appropriate container through a separate waste disposal tube.

2.4.1.11 Sample conveyor of the SPS system

Warning

It is forbidden to touch the sample conveyor of the SPS system during the operation of the analyzer. Otherwise, injury or damage to the analyzer may occur.

Components


Fig. 2.4-15. Schematic drawing of the sample conveyor of the SPS system

(1) Anterior transporter of the chemiluminescence immunoassay module

(2) Posterior transporter of the chemiluminescence immunoassay module

(3) SPS system transfer device

(4) Conveyor support of the chemiluminescence immunoassay module

Function

The conveyor of the SPS module can be divided by function into a conveyor for standard samples, a multifunctional STAT conveyor and a conveyor for returning tripods.

The sample conveyor of the SPS system has the following functions.

- (1) The function of moving the tripod for samples. The sample conveyor of the SPS system can move the sample stand to the appropriate sampling position by the sampler holder.
- (2) Tripod return function for samples. The sample conveyor and the SPS system replacement device allow you to move the tripod with loaded samples to the return area, and then using the transport device back to the sample area. Thanks to this function, the tripod with samples can be returned to the sample area.
- (3) The function of storing a tripod for samples in the buffer zone. There is a buffer zone on the front conveyor of the chemiluminescence immunoassay module, which allows storing a tripod with samples.

2.4.2 SAMPLE AREA

2.4.2.1 Sample area

Note:

The sample area is at room temperature, without the possibility of cooling. For special sample requirements, testing must be carried out in accordance with the reagent instructions in order to avoid affecting the test results.

Components


Fig. 2.4-16. Schematic drawing of the sample area

(1) Sample area framework

(2) Transportation device

(3) Positions for tripods with samples

Function

The sample area is designed to store samples in a sample stand.

With the help of a transportation device, a tripod for samples with test samples is moved to a certain place on the rails.

2.4.2.2 Tripod for samples

Components


Fig. 2.4-17. Schematic drawing of a tripod for samples

(1) Handle

(2) Test tube holder

(3) Identification label of the sample stand

(4) Test tube

(5) Label indicating the type of sample stand

(6) A hole in the base for fixing sample stand

Container

Each sample stand has 10 positions for test tubes. The supported tube types are listed below.

Standard cuvette: $\phi 12 \times 37$ mm, 0.5 ml.

Original sample tube and plastic test tubes: diameter 10-13 mm, height 60-120 mm.

Type

Type	Color of the tube holder	Appointment
Standard tripod	Gray	It is intended for carrying out analyses of standard samples and outpatient samples
Tripod for STAT samples	Red	Designed for urgent sample analysis (STAT)
Tripod for LC inspection	Purple	Designed for LC checks
Tripod for QC	White	Designed for the analysis of control materials

Requirements for the barcode of test tubes

- Type: Code 128, Code 39, Code 93, Codabar, 2/5 Interleaved.
- Resolution: at least 0.007 inches (0.18 mm).
- Numbers: 4-25.
- Size: width — 8-12 mm, working length — no more than 60 mm, width of margins on the sides: at least 3 mm.

2.5 Control unit

The control unit is a computer with installed software for working with the analyzer. It consists of a monitor, a computer, a keyboard and a mouse.

- Basic computer configuration: CPU frequency — ≥ 3.0 GHz, hard disk capacity — ≥ 500 GB; RAM — ≥ 4 GB; RS-232 interface; USB interface; RJ45 network interface; LCD display; keyboard; mouse.
- Supported operating system: Windows 10.
- Security Software: Basic information security software such as 360 Total Security, 360 Antivirus, Tencent PC Manager and Kingsoft Antivirus are supported.
- Software data transfer interface: wired network interface; data is stored in XML format and in MySQL databases.
- User access: User management function is supported; permission management based on two roles is supported: administrators and regular users.
- Network connection: the connection between the software and the analyzer is carried out over a wired network (TCP/IP), while the analyzer does not need to be connected to a local area network (LAN).

2.6 Labels

2.6.1 Module Labels

2.6.1.1 Chemiluminescence immunoassay module

Label with a symbol indicating protective grounding

Label indicating protective grounding

Location:

- ① One label is on the power supply.
- ② One label is on the housing grounding unit.
- ③ Two labels are on the storage compartment of the cuvette.
- ④ Two labels are on the support of the reagent area.

Label "Emergency stop switch" (Emergency stop button)

The emergency stop button.

The label is located under the emergency stop button.

Emergency stop switch

Label "Cuvette bin" (Container for used cuvettes)

The cuvette container is designed to store used cuvettes.

The label is located on the container for used cuvettes.

Cuvette bin

"Tip bin" label (Container for used tips)

The tip container is designed to store used disposable tips.

The label is located on the container for used tips.

TIP bin

Label "Starter box" (Container for starting reagents)

The starter reagent container is designed to store starter reagents.

The label is located on the panel of the container for starting reagents.

Starter box

Label "S₁"

Indicates the location of the starting reagent 1.

The label is located on the side panel inside the container for starting reagents.

Label "S₂"

Indicates the location of the starting reagent 2.


The label is located on the side panel inside the container for starting reagents.

Connector label for connecting an external liquid system

The connectors are designed to connect an external device for clean water, high and low concentration liquid waste.

The label is located above the connectors for connecting the external liquid system.

Label "Test counter" (Analysis counter)

The device is designed to record the number of analyses.

The label is located in the upper left corner of the case.

Test counter

Label "Cuvette storage switch" (Button of the sliding compartment for storing cuvettes)

The button is designed to control the sliding compartment for storing cuvettes.

The label is located on the front panel of the sliding compartment for storing cuvettes.

Cuvette storage switch

Label "Tip storage switch" (Button of the sliding compartment for storing tips)

The button is designed to control the sliding compartment for storing the tips.

The label is located on the front panel of the retractable tip storage compartment.

Tip storage switch

Label "Tip bin switch" (Container button for used tips)

The button is designed to control the container for storing the tips.

The label is located on the left corner of the case.

TIP bin switch

Label "Lower recovery area" (Lower return compartment)

Indicates the lower return compartment in the cuvette storage compartment.

The label is located in the following places:

- ① The right side of the switch of the lower compartment for return;
- ② There is one sheet metal position each to the left and right of the lower return compartment.

Lower recovery area

Label "Lower storage area" (Lower storage compartment)

Indicates the lower compartment of the cell storage compartment.

The label is located in the following places:

- ① The right side of the switch of the lower storage compartment;
- ② There is one sheet metal position each to the left and right of the lower storage compartment.

Lower storage area

Warning label "Do not touch while operating" (Do not touch while operating)

Warns the user not to touch moving parts while working.

Location:

One label is in the storage area under the storage compartment of the cuvette, one label is in the lower compartment for return.

Do not actuate
during operation

Warning label

This label indicates a safety warning; it is located in those areas of the analyzer that the user may be injured in contact with, including the following areas:

- ① There are reagent areas on the inside of the door (two labels).
- ② There are cuvettes on the top panel of the sliding storage compartment.
- ③ On the top panel of the sliding compartment for storing tips.

Label "Switch off main power before opening cover" (Turn off the main power switch before opening the cover)

This symbol and text remind the user to close the analyzer during its normal operation.

Location:

- ① Next to the front cover lock.
- ② Next to the back cover lock.

Please lock the cover when
the analyzer is running.
Switch off the main power
before opening the cover

Label "Close cover during operation" (The lid must be closed during operation)

This symbol and text remind the user to close the insulating cover of the reaction disk during normal operation of the analyzer.

The label is located in the center of the insulating cover of the reaction disk.

Please make the plate cover be closed when the analyzer is working

The label with a warning about opening the case

This symbol and text warn the user that in order to avoid electric shock, it is forbidden to open the analyzer housing. If maintenance is required, contact our company or an authorized maintenance specialist. The label is located on the back of the case.

To avoid electric shock, do not open the shell. For service, please contact our company or the authorized professional personnel.

The label warning of the danger of biological contamination

Reminds the user of the risk of biological contamination in contact with the relevant areas of the analyzer.


Location:

- ① On the handles of containers for used cuvettes and tips.
- ② Above the liquid waste outlet on the back of the analyzer.


Label "Thick Waste Liquid" (High concentration liquid waste)

This symbol and text remind the user that high-concentration liquid waste is stored in the appropriate tank. It is located on a high-concentration waste tank.

Label "Thin Waste Liquid" (Low concentration liquid waste)

This symbol and text remind the user that low-concentration liquid waste is stored in the appropriate tank. It is located in a tank for low-concentration waste.

Thin Waste Liquid

Label "Washing liquid"

Indicates the container in which the washing liquid is stored.

It is located on the tank for loading the washing liquid and the capacity for fixing the washing liquid.

SYSTEM
LIQUID

Label "Pure Water" (Pure water)

Indicates the container in which clean water is stored.

It is located on a container for clean water.

PURIFIED
WATER

Label "No mixing" (Do not mix)

Reminds users that starting reagent 1 and starting reagent 2 cannot be mixed.

The label is located on the inside of the lid of the container for starting reagents.

**2.6.1.2 Sample area****The label warning of the danger of biological contamination**

Reminds the user of the risk of biological contamination in contact with the relevant areas of the analyzer.

There are reagent areas on the left support.

2.6.2 Other labels

Label	Description
	Manufacturer
	Authorized Representative in the European Community/The European Union
	Production date
	Catalog number
	Medical device for in vitro diagnostics and monitoring
	Serial number
	The warning. Refer to the attached documents


WEEE Label

The presence of the WEEE label on the product means that it should not be recycled as household waste. Proper disposal of this device will help to avoid the effects of hazardous substances on the environment and human health. For more information, contact the manufacturer or distributor from whom you purchased this equipment.

This side up

This symbol indicates that the packaging should be placed vertically during transportation. This symbol is located at the top in the center on each of the four sides of the package.

Keep away from moisture

This symbol indicates that the packaging must be protected from moisture during transportation.

This symbol is located at the top in the center on each of the four sides of the package.

Fragile products

This symbol indicates that the contents of the package are fragile and must be handled carefully.

This symbol is located at the top in the center on each of the four sides of the package.

Don't roll

This symbol indicates that it is forbidden to roll the box during transportation.


This symbol is located at the top in the center on each of the four sides of the package.

Stacking is prohibited

This symbol indicates that the transport package should not fit in more than one tier.

This symbol is located at the top in the center on each of the four sides of the package.

2.7 Scope of delivery

See the packing list.

2.8 Required materials (not included in the package)

Title	Catalog numbers	Note	NO. RU
In vitro Reagent Kit for automatic analyzers MAGLUMI®	individually	/	individually
Reaction Cups for automatic analyzers MAGLUMI®	130105000101	/	/
Tip for automatic analyzers MAGLUMI®	130208000201	/	/
Waste Bag for automatic analyzers MAGLUMI®	130210000201	Suitable for a container for used cuvettes	RZN 2021/14988
	130210000101	Suitable for a container for used disposable tips	
Wash Concentrate for automatic analyzers MAGLUMI®	130299005M	With manual preparation	RZN 2021/14863
	130299035M	In case of automotive training	
System Tubing Cleaning Solution for MAGLUMI® automatic analyzers	130299007M	/	RZN 2021/14881
Starter Reagent Kit (Starter 1, Starter 2) for automatic analyzers MAGLUMI®	130299027M	/	RZN 2021/14858
Luminescence measurement control (Light Check) for MAGLUMI® automatic analyzers	130299006M	/	RZN 2021/14861

Chapter 3 ANALYZER INSTALLATION

3.1 Storage and transportation requirements

3.1.1 Storage requirements

The packaged Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories must be stored in well-ventilated rooms in which there is no corrosive gas, at temperatures from -20 to 55 °C, relative humidity no more than 93% and atmospheric pressure from 50.0 to 106.0 kPa.

3.1.2 Transportation requirements

The Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories must be transported in a package at an ambient temperature from -20 °C to 55 °C and a relative humidity level of no more than 93%. It is necessary to protect the analyzer from shocks, moisture and sunlight.

3.2 Installation Requirements

3.2.2 Requirements for spatial spacing

To ensure personnel access to the analyzer during operation, maintenance and repair, its installation must be carried out in accordance with the following requirements:

- The distance between the wall and the left side of the analyzer should be at least 80 cm.
- The distance between the wall and the right side of the analyzer should be at least 150 cm.
- The distance between the back panel of the analyzer and the wall should be at least 80 cm.
- The distance between the front panel of the analyzer and other devices should be at least 130 cm.
- There should be enough space to install a liquid waste drain device.
- There should be enough space near the outlet to properly connect and disconnect the power cord. The analyzer should be installed in such a way that it can be easily disconnected from the power supply.

3.2.3 Environmental requirements

- Ambient operating temperature: 10-30 °C
- Relative humidity: ≤70 %
- The room should be well ventilated, free of dust, mechanical vibrations, loud noise and interference from the power supply network.

- Do not place the analyzer near an electric brush motor, a source of flashing fluorescent light, and electrical devices that turn on and off frequently.
- Avoid direct sunlight and do not place the equipment near heat sources and air flows.
- Avoid direct sunlight and do not place the equipment near heat sources and air flows. The analyzer should be installed on a flat surface that can withstand a load of at least 1049 kg.

3.2.4 Power supply requirements

- This analyzer requires a 100-240 AC power supply At a frequency of 50/60 Hz.
- Rated power of the analyzer: 2200 V·A
- Overheating protection device: 2455R-70820258 (Technical specifications: IN THE OPEN STATE — $53^{\circ}\pm 3^{\circ}$; IN THE CLOSED STATE — $45^{\circ}\pm 4^{\circ}/-5^{\circ}$, DC 12V, 15A), 9700K51-1638 (technical specifications: IN THE OPEN STATE — 75°C ($\pm 5^{\circ}\text{C}$), IN THE CLOSED STATE — 50°C ($\pm 15^{\circ}\text{C}$), 230 V AC current, 13 A).
- For the smooth operation of the analyzer, it is necessary that the value of the backup power of the distribution unit is more than 20%, and the power of the uninterruptible power supply is at least 5 kVA.
- To ensure safety, the analyzer must be connected to a properly grounded power source.
- For safety reasons, an overcurrent protection device with a rated current of 15 A must be installed between the building's power supply and the device's power cord.

Warning

Failure to comply with the above environmental and power supply requirements may adversely affect the operation of the analyzer, cause damage to it, or even lead to personal injury.

3.2.5 Water supply requirements

- The conductivity of pure water is less than 1 $\mu\text{S}/\text{cm}$, the pressure of purified water should be from 0.1 MPa to 1 MPa, and the flow rate of pure water should be greater than 10 l/h.

Warning:

The water supply for the analyzer must meet the above requirements for clean water, otherwise the results may be inaccurate or the analyzer may not work correctly.

3.2.6 Drainage requirements

- **Wastewater is a variety of liquid waste generated during the operation of the device, which is divided into concentrated and diluted waste:**

High concentration liquid waste is higher concentration waste, which mainly comes from the sample arm, measuring chamber and cleaning station, the degree of contamination is high;

Low-concentration liquid waste is all liquids except concentrated waste, including reagents, washing liquids, condensate, with the exception of samples, the degree of contamination is lower.

- Waste products can be stored in a special external waste barrel or can be processed by the user in accordance with local requirements for the treatment of medical waste.

**Warning:**

The discharge of waste liquids must be carried out in accordance with local regulations for the disposal of medical waste.

3.3 Unpacking and checking

3.3.2 Unpacking procedure

After receiving the analyzer, carefully inspect the external surfaces for damage. In case of damage, please contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or to a local distributor. After you have verified that there is no damage, follow the steps below to open the package.

- The package must be placed vertically so that the arrow on the label located on the package points upwards.
- Open the package and check if all the components listed on the packing list are present. If any component is missing, please contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or to a local distributor.
- Carefully inspect the external surfaces of the analyzer for damage. In case of damage, please contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or to a local distributor.

3.3.3 Movement and location

- When moving short distances, the analyzer can be pushed (the device is equipped with built-in wheels). Remove the feet before moving.
- When moving the analyzer, keep the analyzer upright, do not tilt or place the device on its side. Move the analyzer carefully and slowly.
- Do not allow vibration during handling. All detected faults must be checked and eliminated before moving and operating the device.
- The four legs of the analyzer should be lowered to the floor level so that they securely fix the device after it is placed in the required position. The analyzer should be securely fixed, should not wobble, and the length of the legs should be adjusted so that the analyzer is level.

3.4 Analyzer Installation

Installation of the analyzer should be performed by our company's specialists. It is not recommended that the analyzer be disassembled and installed for purposes other than standard system maintenance.

3.4.2 Connection between modules

- The sample area and chemiluminescence immunoassay module are secured by two guide studs and four hexagonal head cap screws (M12×55).

- Connect the power bus and CAN bus between the modules.

3.4.3 Electrical connection

- Connect the monitor, keyboard, and mouse to the appropriate connectors on the computer.
- Connect the monitor cable to the USB port on the back of the computer.
- Connect the power cables of the system unit and monitor to the appropriate ports.
- Connect the analyzer's power cable to a connector that meets the specified requirements.
- Connect the appropriate connectors of the three network ports on the router separately to the computer network port, the control panel network port on the back of the analyzer, and the internal network port of the sample area.

3.4.4 Connecting washing liquid containers

By default, the analyzer uses automatic flush fluid preparation mode, if manual flush fluid preparation mode is used, open the analyzer enclosure door, remove the flush fluid tray, remove the clean water container, Connect the fluid suction component and connecting tube to the outlet tube of the bottom of the clean water container.

3.4.5 Connection of liquid waste tanks

If using an external high concentration liquid waste tank, connect the high concentration liquid waste outlet marked (THICK_WASTE) on the back of the analyzer and the high concentration liquid waste tank inlet using the appropriate connectors. Then connect the High Concentration Liquid Waste Sensor to the High Concentration Liquid Waste Sensor connector.

If using an external low-concentration liquid waste tank, connect the low-concentration liquid waste outlet marked (THIN_WASTE) on the rear panel of the analyzer and the low-concentration liquid waste tank inlet using the appropriate connectors. Then connect the Low Concentration Liquid Waste Sensor to the Low Concentration Liquid Waste Sensor connector.

In the case of custom waste discharge tubes, use the tubes to connect the outlet ports for high and low concentration liquid waste to the discharge tubes separately.

3.5 Starting and checking the system status

After installation, make sure that the following elements are properly installed and securely connected: power cables and network cables of the analyzer and computer, flushing liquid tubes, high and low concentration liquid waste systems. Turn on the power and start the system.

1. Turn on the power of the Fully-auto Chemiluminescence Immunoassay Analyzer MAGLUMI X8 with accessories and turn the switches "UNIT" (Device) of the immunochemiluminescent analysis module and "RACK" (Tripod) of the sample processing system located on the control

- panel on the back of the analyzer to the "on" position.
2. After logging in to Windows, launch and log in to the software.
3. After logging into the software and running the user version, the analyzer will enter the initialization state and initialization will be performed for each component.
4. Load consumables such as cuvettes, disposable tips, flushing liquid and starting reagents.
5. Ensure that the status of each component is displayed in the software, temperature, humidity and electrical voltage are within normal ranges, check that waste and solid and liquid consumables are adequate and normal.
6. Rinse the sample dispenser holder, reagent dispenser holder, flushing unit and chamber tubes completely.
7. Follow the system verification procedure. Perform 10 BGW checks on all chemiluminescence immunoassay modules, 8 LC checks on the sampler and 5 LC checks on the left and right reagent dispensers (instructions and information about the expected results are given in the section 4.2.4 "Checking the analyzer status"). Upon receipt of the system verification results that meet the requirements, the verification is completed and the installation of the analyzer is confirmed.

Chapter 4 BASIC OPERATING PROCEDURES

This section describes the basic operating procedures of the system. After studying this section, the user will be able to perform standard operating procedures using the system. For detailed instructions on how to use each interface and function of the user software, see chapter 5 of this manual.

4.1 Standard operating procedures

Table. 4.1-1. Standard operating procedures

The procedure of actions	Procedure
1. Pre-launch checks (1) Checking the power supply (2) Checking the analyzer	Check the operability of the power supply, the analyzer's water supply meets the requirement Check the appearance and status of the analyzer according to the instructions
2. Connecting the power supply	Turn on the power supply and turn the switches "UNIT" (Device) of the chemiluminescence immunoassay module and "RACK" (Tripod) of the sample processing system to the on position
3. Logging in to the user software	Enter the operator's username and password in the software login dialog box.
4. Replacement of consumables	Make sure that the volume of the flushing liquid, starting reagent 1, starting reagent 2, the number of cuvettes and disposable tips is sufficient for the current analysis. If there are not enough supplies, replenish their quantity in a timely manner.
5. Checking the analyzer status (1) Status check (2) Confirmation of the results of the system check	Wait at least 30 minutes after turning on the analyzer to make sure that it is in a stable condition. Click on the icon of the desired tool in the main menu to go to the status interface of this tool. Check the environment, waste, consumables, tools and handling of substances. Before starting the analysis, make sure that all states are displayed in green. Follow the system verification procedures and make sure that the measurement values of the LC and BGW tests are within the acceptable range
6. Checking the conditions for performing the analysis (1) Checking the analysis parameters (2) Profile verification	Check the analysis parameters. Configure the profile according to the current conditions

7. Loading of reagents (1) Reagent loading (2) Checking the remaining reagent volume (3) Checking the mixing time of magnetic microparticles	Load the reagent properly. Make sure that the number of remaining reagent tests is sufficient for the analysis. Make sure that the magnetic microparticles of the reagent have been stirred for at least 30 minutes
8. Register the calibration analysis. Register the QC analysis.	Make sure that the reagent has a valid calibration record. Perform calibration for reagents that do not have valid calibration data. Do the following The QC of the analysis and make sure that the position of the QC liquid registered for the QC matches the liquid on the appropriate sample stand.
9. Uploading samples	Upload a sample for testing in the sample area.
10. Sample registration	Register a standard sample. Register an urgent test (STAT). Register the diluted sample
11. Starting the analysis	After you confirm that the registration list is correct, click <Start> on the action bar to start the analysis.
12. Analysis of added samples	Register all samples added during the analysis. (if the sample is a STAT sample, it must be registered accordingly; if the sample is diluted, it must be registered accordingly). Click <Start> on the Quick Access toolbar to start performing the analysis.
13. Track the performance of the analysis. (1) Status Tracking (2) Suspension of transportation (3) Stopping machine operation	Change module status in real time, including environmental and consumable conditions To pause transport of the sample rack during analysis, press <Suspend Transport>. To stop the current analyzer action while the analysis is in progress, click <Stop Analyzer>.
14. Working with results	Query the analysis results, re-calculate and re-measure, delete or print the results
15. Until the procedure is completed	Remove the reagent kits from the reagent area and the samples from the sample area. Empty the waste, used cuvette, and tip containers.
16. Upon completion (1) Turn off the power to the analyzer (2) Place the analyzer in sleep mode	Exit the software and turn off power to the analyzer and computer. If you do not want to turn off the analyzer and computer, put them into sleep mode ([Home Screen] > <Auto Hibernate >).

4.2 Preparation for analysis

Follow the steps below to prepare for the standard assay.

4.2.1 Pre-launch checks

To ensure proper system operation, the user should check the following aspects in advance.

1. Check the power supply to ensure normal power supply.
2. Check the network data cable between the analyzer and computer and make sure they are properly and securely connected.
3. Make sure flushing fluid and starting reagent reservoirs 1 and 2 are properly connected and contain enough fluid.
4. Make sure that the number of cuvettes and disposable tips is sufficient.
5. Make sure the liquid waste tank is empty.
6. Ensure that used cuvette and tip containers and spent cuvette and tip assemblies are empty.

Warning

Be aware of the risk of infection from contact with biological liquids. To avoid contamination, always wear gloves and special clothing when performing the tests described; if necessary, wear protective eyewear.

4.2.2 Turning on and logging into the software

1. Turn on the power supply and turn the switches "UNIT" (Device) of the chemiluminescence immunoassay module and "RACK" (Tripod) of the sample processing system to the on position
2. After logging into Windows, double-click the user software shortcut on your desktop to launch the software. After starting the software, a login dialog box will appear on the screen. Enter the user name and password and click <Login> to enter the user software interface.
3. After the user logs into the software, the analyzer will enter the initialization state and initialization will be performed for each component. You can run tests after all modules have been placed in standby mode.

Note.

1. The system administrator account username is "snibe" and the original password is "snibe".
2. To access the [Home Screen] page, the user software must be connected to the module, otherwise the message "Analyzer Search" will be displayed after logging into the software and attempting to access the main menu page.

4.2.3 Replacement of consumables

If the analyzer does not have enough cuvettes, wash fluid, disposable tips, or starting reagents, or if

any of the consumables are missing, an alarm will sound. The analyzer will return to normal after the missing consumables have been loaded.

4.2.3.1 Loading cuvettes and disposable tips

● Loading cuvettes

New and used cuvette blocks are stored in the sliding cuvette storage compartment. The sliding compartment is located on the left side of the analyzer housing, and each movement can be controlled by the corresponding sliding compartment button, lower storage compartment button, and lower return compartment button.

Button	Indicator light color	Description
Sliding compartment button	Green	The sliding cuvette storage compartment is in standby mode. When the button is pressed, the sliding compartment will open or close by itself
	Flashing green	The sliding cuvette storage compartment is in the startup state. To stop movement of a sliding compartment immediately, click the button for that compartment. Click the sliding compartment button again to continue moving
	Orange	The sliding cuvette storage compartment is in the "occupied" state and the button is inaccessible
	Flashing orange	Cuvette sliding storage compartment error, button not available. In this case, push (or pull) the sliding compartment so that it starts to move inward (or outward)
	Doesn't light up	The analyzer has not been initialized, so the sliding cuvette storage compartment can be opened or closed manually
Button bottom storage compartment	Green	The lower compartment of the cuvette storage compartment is in standby mode and ready for use. When the button is pressed, the height of the lower cuvette storage compartment will be adjusted to the most suitable height to accommodate the cuvette blocks
	Orange	The lower cuvette storage compartment is in the "occupied" state and the lower storage compartment button is inaccessible
	Doesn't light up	The analyzer has not been initialized, so the lift mechanism of the lower storage compartment can be moved by pulling it manually
Bottom compartment button	Green	The lower return compartment in the cuvette storage compartment is on standby and ready for

for return		use. When the button is pressed, the height of the lower return compartment will be adjusted to the most suitable height for retrieving cuvette blocks
	Orange	The lower return compartment in the cuvette storage compartment is in the "occupied" state and the lower return compartment button is not available
	Doesn't light up	The analyzer has not been initialized, so the frame lift mechanism of the lower return compartment can be moved by pulling it manually

The procedure of actions%

1. Log in to the user software. Once the analyzer has completed initialization, check the software for the number of new cuvettes and used cuvettes available.
2. Remove the cuvettes from the packaging and prepare the number of new cuvette blocks required for analysis.
3. Wait until the drawer button light is green, and then press the drawer button. The sliding compartment will then slide out of the analyzer housing and the button indicator will flash green.
4. Wait until the sliding compartment button and lower storage compartment button lights turn green, then load a new cuvette assembly.
5. Grasp one or more cuvette blocks by the handle with both hands and slowly place the blocks into the bottom compartment of the sliding compartment (see Figure 4.2-1). Press the button of the lower storage compartment and its height will be adjusted to the loading level. Wait until the bottom storage compartment button indicator light turns green before installing the following cuvette blocks.
6. After the lower return compartment button indicator light turns green, grasp one or more spent cuvette blocks by the handle with both hands and slowly remove the blocks from the lower return compartment of the pull-out compartment. Press the lower compartment button to return, then its height will adjust to the load level. Wait until the Lower Return Compartment Button LED turns green before removing the following cuvette blocks.
7. After replacing the cuvette assembly, click the sliding compartment button. The sliding compartment button indicator light will start flashing green. Wait for the sliding compartment to move inside the analyzer to the proper position and for the button indicator to turn green. After that, the download procedure is considered to be completed.


Fig. 4.2-1. Installation of a new cuvette assembly

Note.

1. Pay attention to the number of remaining cuvettes before and during analysis.
2. If the analyzer has not been used for one week, empty the cuvette storage compartment before use to ensure operation with clean cuvettes.
3. The cuvette blocks and cuvette storage compartment are designed to avoid improper installation. Install these components in the correct direction.
4. After installing a new cuvette block, be sure to remove the spent cuvette blocks or normal analyzer operation may be impaired.

● **Loading disposable tips**

The sliding tip storage compartment for loading disposable tips is located on the left side of the analyzer housing. The various movements of the sliding compartment can be controlled using the buttons located on the compartment.

Button	Indicator light color	Description
Sliding compartment button	Green	The sliding tip storage compartment is on standby and ready for use. When the button is pressed, the sliding compartment will open or close by itself
	Flashing green	The sliding tip storage compartment is in the startup state. To stop movement of a sliding compartment immediately, click the button for that compartment. Click the sliding compartment button again to continue moving
	Orange	The sliding tip storage compartment is in the "occupied" state, the button is inaccessible

	Flashing orange	Tip storage sliding compartment error, button not available. In this case, push (or pull) the sliding compartment so that it starts to move inward (or outward)
	Doesn't light up	The analyzer has not been initialized, so the sliding tip storage compartment can be opened or closed manually

The procedure of actions%

1. Log in to the user software. After the analyzer initialization is complete, check the software for the number of disposable tips available.
2. Remove the disposable tip assembly and prepare the number of disposable tips required for the assay.
3. Wait until the drawer button light is green, and then press the drawer button. The sliding compartment will then slide out of the analyzer housing and the button indicator will flash green.
4. Wait for the sliding compartment button indicator light turns green, then you can begin loading the new tip assembly. And remove the tip assembly with the used disposable tips
5. Grasp both sides of the new tip assembly and place in a free position in the sliding compartment (see Figure 4.2-2).


Fig. 4.2-2. Loading new disposable tips

6. After loading the tip assembly, click the sliding compartment button. The sliding compartment button indicator light will start flashing green. Wait for the sliding compartment to move inside the analyzer to the proper position and for the button indicator to turn green. After that, the download procedure is considered to be completed.

**Note.**

1. Pay attention to the number of disposable tips remaining before and during analysis.
2. If the analyzer has not been used for one week, empty the tip storage compartment before use to ensure operation with clean disposable tips.

4.2.3.2 Adding washing liquid and starting reagent

- **Adding washing liquid**

The container with the wash fluid solution is placed on a tray in the center of the analyzer cabinet. The remaining fluid volume can be monitored on the home screen of the respective analyzer.

The procedure of actions%

1. Open the front door of the analyzer enclosure, slide out the flushing fluid tray, and remove the suction assembly from the upper flushing fluid reservoir. Then remove the empty flushing fluid reservoir.
2. Place a new flushing fluid reservoir on the front of the tray, unscrew the top lid of the flushing fluid reservoir, and screw the lid on and then insert the flushing fluid reservoir suction assembly. Insert the flushing fluid tray back into the analyzer.

Note.

- 1) When using the manual mode, the flushing liquid must be prepared in advance. Ensure that the flushing fluid reservoir is free of contaminants prior to preparation.
- 2) Pay attention to the volume of remaining wash fluid before and during analysis.

- **Adding the starting reagent**

The procedure of actions%

1. Slide out the starter reagent container and lift the starter reagent tube support to remove the empty starter reagent container.
2. Place the new starter reagent container in the proper position, lower the starter reagent tube support, and slide the starter reagent container into the analyzer housing.

Note.

1. If the chamber tubes are empty before flushing, flush the chamber 50 times before performing a BGW check.
2. Pay attention to the remaining volume of starting reagent before and during the analysis.
3. The containers for the starting reagent 1 and 2 as well as the tubes are labeled.

Do not mix up the components during installation!

4.2.4 Checking the analyzer status

4.2.4.1 Status check

Tap the appropriate analyzer icon on the home screen or in the main menu to go to the home screen of that analyzer. Check the environment, waste, consumables, tools and handling of the substance. In case of correct status, the icons are displayed in green color.

After the analyzer initialization is complete, wait 30 minutes to ensure that the environmental parameters are stable.

If the waste status is outside the normal range, follow the steps described in Section 4.6.5, Waste Disposal, of this document.

If the status of the consumables is out of the normal range, follow the steps outlined in Section 4.2.3, "Replacing Consumables" of this document.

If the status of the tool is displayed as unacceptable, contact our service technicians.

If the transaction status is displayed as invalid, perform the appropriate calibration, system check, or maintenance procedure.

Note.

After all analyzer components have completed initialization, wait for all status lights to turn green before proceeding with the analysis. If the analyzer does not return to normal operation for an extended period of time, please contact our service technicians.

4.2.4.2 Confirmation of system test results

The system check consists of a background flush (BGW), light check (LC) for the sampler and reagent dispenser (left and right reagent dispenser).

Before performing a system check, place the light test solution in the appropriate positions in the sample area and the right and left reagent areas. See 5.6.2 "Loading liquid to check the light signal" for details.

Press the <Utility> button on the main menu to go to the [Utility] screen. Click <Add> to open the [Registration of data for system verification] screen, enter the sample arm, sampler, number of reagent dispenser flushes, flush assembly, and chamber, and the number of BGW checks and LC checks for the sampler and reagent dispenser (left and right reagent dispenser). Click <Registration> and then <OK> to start the system check. When finished, check the results of the system check.

Perform 10 BGW checks, 8 LC sampler checks, and 5 LC reagent dispenser checks. The system test is considered passed when results are obtained that meet the following requirements.

● BGW	100–1200	coefficient of variation $\leq 10\%$
● LC	450000–650000	coefficient of variation $\leq 3\%$

Note.

1. Analyses may be performed only if the results and coefficients of variation of the BGW and LC checks are satisfactory. It is not recommended to perform sample analysis if the test results are outside the acceptable range (deviating values are highlighted in red in the system test results list). Otherwise, it may adversely affect the validity of the assay results.
2. The standard LC check should be performed once a week. The LC test should also be performed when analyzer parts are replaced. If the LC test is not performed, fluid loading is not required to test the light signal.
3. BGW flushing should be performed every day prior to analysis to ensure accurate results.

4.2.5 Checking the conditions for performing the analysis

Ensure that the parameters and reference range values of the assay are set correctly before performing the assay. Ensure that the assay reagent is calibrated and the calibration result is approved and not expired. Customize profile, dilution, and material data for QC.

For detailed information on setting the assay, dilution, and profile parameters, see 5.7.1 "Setting Parameters".

For detailed information about setting up materials for QC, see the 5.5.1 "Registration of QC" section.

4.2.6 Loading of reagents

Scan the RFID tag of the reagent by positioning it between the two reagent areas, then place the reagent in the reagent area. Check the following information on the [Reagent] screen of the software before calibrating, analyzing control materials, and analyzing a sample using reagent.

- Ensure that the reagent has not expired.
- The amount of remaining reagents should be sufficient to perform the assay.
- Magnetic microparticle reagents should be stirred for 30 minutes.

4.2.6.1 Design of the reagent kit

The reagent kit consists of the following components.


Fig. 4.2-3. Design of the reagent kit

1. Sealing film: serves to seal the reagents. The analyzer reagent dispenser is inserted into the reagent kit through a sealing film for reagent withdrawal.
2. Handle: used to load and remove the reagent kit.
3. RFID tag: contains information about the reagent kit.
4. Identification label: contains information about the components, name, expiration date, etc. of the reagent kit.
5. Fixing clip: designed to fix the reagent kit when inserted into the reagent area.
6. Identification blocker: Designed to block the optocoupler when installing a reagent kit in the reagent area for its identification.
7. Mechanism for stirring magnetic microparticles: together with the shaker of the reagent area, the magnetic microparticles are stirred.

4.2.6.2 Loading the reagent kit

1. Preparation of the reagent kit

Ensure that the reagent kit is stored in the correct position (sealing film should be facing up). Do not invert or shake the reagent kit.

- 1) Remove the reagent kit package from the refrigerator.
- 2) Remove the reagent kit from its packaging.
- 3) Ensure that the reagent kit is not leaking. If leakage is present, contact Shenzhen New Industries Biomedical Engineering Co., Ltd. or distributor in a timely manner.
- 4) Gently rotate the magnetic microparticle stirring mechanism of the reagent kit to ensure that nothing is obstructing the movement of the mechanism.

- 5) Carefully remove the aluminum film from the reagent kit.
- 6) Check the fluid in each component of the reagent kit for air bubbles. If air bubbles are present, use a pipette to remove them. Ensure that bubbles are completely absent before use.
- 7) Remove fluid from the surface of the reagent kit to avoid cross-contamination.

Note.

- 1) Remove the protective foil from the reagent kit before using a new reagent.
 - 2) Reagents should be prepared, used and stored strictly in accordance with the reagent-specific instructions for use. Do not allow air bubbles to form in the reagents. This may affect the accuracy of reagent collection and analytical results.
 - 3) Reagents from different kits should not be mixed as this may affect the validity of the assay results.
-

2. Loading the reagent kit

Ensure proper operation of the analyzer and software before inserting the reagent kit into the reagent area.

- 1) Open the door of the reagent area.
- 2) Remove unused reagent kits according to the indicator lights of the various channels in the reagent area.
- 3) Holding the reagent kit by the handle, bring the reagent kit with the RFID tag toward the readout area (for approximately 2 seconds). If scanning is successful, a beep will sound.
- 4) Insert the reagent kit along the empty reagent well. When the reagent kit is installed correctly, the analyzer will beep.
- 5) The reagent data will be read by the software. The corresponding item will be highlighted on the [Reagent] screen, where the name of the reagent will also be displayed.
- 6) If an error occurs while reading data, repeat the above steps.
- 7) When you have finished reading the reagent information, close the reagent area door.

The status of the reagent kit is indicated by color indicators in the reagent area.

Indicator light color	Reagent kit status
Doesn't light up	Standby mode, reagent kit can be loaded
Lights up green	Standby mode. The reagent kit is installed in a stand
Flashing green	The dispensing of the reagent kit is complete and the reagent rack can be removed.
Lights up yellow	Startup status. The reagent kit is used
Flashing yellow	Error. For example, the reagent kit was removed during use

3. Checking the remaining volume of reagent

After the RFID tag has been scanned into the reagent area and the reagent kit has been properly installed, you should check the number of remaining tests in the reagent area and verify that the remaining tests are sufficient for this assay. If the reagent volume is insufficient to complete the current sample analysis, install a new reagent kit in the reagent area.

4. Checking the mixing time of magnetic reagent microparticles

After the reagent set is installed, the [Reagent] dialog box will display the reagent set name, the number of remaining tests for the reagent, the calibration status, and the magnetic microparticle agitation time in the corresponding reagent rack area. See the Reagent Information section for details.5.4.1

4.2.7 Uploading samples

4.2.7.1 Sample preparation

Ensure that the following conditions are met before analyzing the sample.

- 1) Do not use samples with hemolysis or lipemia or samples contaminated with microorganisms for analysis.
- 2) Ensure that there are no air bubbles in the sample before analysis.

Verify that the sample is appropriate for the conditions, and then insert the tube into the sample rack. The procedure is described below.

- 1) The test tube shall meet the requirements described in Section 2.4.2.2 Sample Rack.
- 2) Carefully insert the test tube into the sample rack.
- 3) If there is a barcode on the tube, it must be turned with the barcode against the open side of the rack so that the barcode scanner can read the required data.

Note.

1. To prevent damage to the sampler holder, make sure the tube lids are open before adding the sample.
 2. To ensure reliable results, sample collection and handling must comply with all listed requirements and performance characteristics to avoid air bubbles and blood clots.
 3. Do not rotate the tubes with the labeled barcode in the sample rack as this may damage the barcode.
-

Warning

1. Wear gloves when handling samples to avoid infection.
 2. If the sample volume is insufficient or no sample is available, the analysis
-

result displays a concentration value labeled "L", "S", or other. If the analyzer displays the message "The sampler pipette cannot draw a sample", the analysis must be run again.

3. After placing the test tube rack in the sample area, do not change the position of the tubes in the rack to avoid confusion of test results.

Installing the Light Check Liquid (LC) supplied by our company into the rack is the same as installing the sample tube, but the barcode on the tube may not be turned toward the barcode scanner when installing the Light Check Liquid. The LC test is only performed when the analyzer recognizes the LC test tripod. If the test tube is not properly installed in the rack, remove it and reinstall it properly.

4.2.7.2 Loading the sample rack

Follow the steps below to load the sample rack.

- 1) Open the door of the sample area.
- 2) Remove unused sample racks according to the color of each rack indicator in the sample area.
- 3) Select an empty position for the sample tripod and gently slide the tripod along the block so that it is fully in position.

Indicator light color	Status of the sample rack
Doesn't light up	Often occurs when sample racks are not loaded. If the sample rack is already loaded, this indicates a failure of the sensor or other equipment.
Lights up green	Waiting
Flashing green	The analyzer has finished loading samples, and the sample rack can be removed
Lights up yellow	The sample rack is loaded
Flashing yellow	Indicates an error: the current position for the sample rack is occupied or the current rack has been removed.

Note.

Ensure that the analyzer and software are working properly before placing the sample rack in the sample area. Otherwise, the sample rack may not be recognized and subsequent operations will be invalid.

4.3 Analyzing

Once the preparatory procedures have been completed, the analysis can begin.

4.3.1 Registration of calibration analysis

From the main menu, click the <Reagent> button to go to the [Reagent] screen. Select the reagent to calibrate and press <Reg. cal. >. Confirm the assay for which you want to calibrate in the dialog box and press <OK> to start the calibration.


Fig. 4.3-1. Reagent Screen

When the calibration is complete, double-click the appropriate reagent to be analyzed or select it and click <Calibration Curve> to go to the [Reagent Information] screen. Check the curve details and accept or reject the new curve.

To approve the new calibration data, select the "Confirm" option and click the <Confirm> button. To reject the new calibration data, select the "Reject" option and click the <Confirm> button. You can re-register the reagent and perform calibration using the [Reagent] screen.


Fig. 4.3-2. Detailed information about reagent calibration

4.3.2 Registration of the QC analysis

Click **<QC>** in the main menu. On the **[QC]** screen, click on the **<Reg. of data for QC>** tab to go to the **[Reg. of data for QC]** screen. After correctly loading the tripod with the materials for the QC, select the materials QC and analyses from the list, and then specify modules and positions for tripods. Click **<Registration>** to complete the registration of the QC tests.


Fig. 4.3-3. Registration of data for QC

4.3.3 Sample registration

Click the **<Sample>** button in the main menu to go to the **[Sample]** screen.

Load the samples for analysis into the sample stand. If there is a barcode on the test tube, turn the test tube with the barcode to the open side of the tripod. Install the sample tripod on the sample platform properly, click **<Scan>** and select the desired tripod to scan. Then click **<Scan>** and the analyzer will scan and identify the barcodes of all tripods in the selected sample area. Upon completion of the scan, the scanned sample tripod numbers and scanned barcode numbers will be displayed in the software.

Select the rack and position for the sample data to be changed, enter the sample ID (no or optional ID if barcode is present), select the desired assay in the assay selection area, or use the profile function to select multiple assays and register samples.

When registering an urgent sample (STAT), install a tripod with a STAT sample immediately in the sample area and scan the barcode. Registration of such samples is performed in the same way as registration of standard samples, however, STAT tripods are given priority when performing analysis.

If the sample needs to be diluted, select the analysis that requires dilution. Then click **<Dilution>**, select the dilution ratio and click **<OK>**.

To use the LIS function, set the settings correctly and connect the LIS server, scan or enter the sample ID. The analyzer will automatically extract the analysis data from the LIS server. Another way: Tap **<Use for FOXES>** to select a sample on the **[Use for FOXES]** screen, and then tap **<Use for FOXES>** to get the sample analysis data.

When registration is complete, you can click **<Worklist>** to view the status of all registered tests, delete or change the number of repetitions, and export the data. You can also click the **<Overview of**

Information> button to view the modified analysis data for each sample, the number of registrations for each analysis, and the number of remaining tests for the reagent.


Fig. 4.3-4. Sample

4.3.4 Starting the analysis

When the sample registration is complete, click the **<Start>** button on the action bar to send the appropriate command; after that, the analyzer will start the analysis procedure.

To use the LIS feature, you can press the **<Start>** button after installing the tripod. The analyzer will automatically perform a number of operations, including scanning the sample tripod, scanning the barcode of the sample, loading the analyses and running the analyses.

Note.

- 1) Do not open the analyzer cover during operation.
- 2) Do not touch the sampler holder, reagent probe, reaction disk or other movable devices, and do not prevent them from moving while the analyzer is operating.

4.3.5 Analysis of added samples

1. Additional analysis of samples with recorded data

If you need to add analyses during testing, you will need to register the sample before adding it. In the sample area of the **[Sample]** screen, select a tripod for samples, in the registration area, select the sample for which you want to add analyzes, and then register the sample. For more information about the registration methods, see in the section 4.3.3 "Sample registration".

2. Adding sample analyses without recorded data

If it is necessary to add a sample during the analysis, then after installing the added samples in the sample stand and properly loading the tripod in the sample area, the added samples should be registered (see section 4.3.3 "Sample registration").

If there are not enough reagents in the reagent area to perform additional analyses, If the reagents of additional analyses are not placed in advance in the reagent area, load the appropriate reagents into this area. If reagents are loaded for additional analyses after registration, click the **<Start>** button on the action bar to complete the continuous loading procedure.

4.4 Tracking the execution of the analysis

During the analysis, you can monitor various analyzer states on the home screen in real time (including the status of consumables, temperature and voltage readings). If you need to pause the output of the sample tripods in the sample area during the analysis, click the button **<Suspend Transport>**. If a serious problem occurs, click the **<Stop Analyzer>** button to stop the analysis. For more information about the states, see the 5.9 "Description of the functions of the [Status] screen" section.


4.5 Working with results

Click the **<Result>** button in the main menu to go to the **[Result]** screen. This screen displays the results of all sample analyses performed during the day.

To check the detailed information about the analysis results, double-click on the result record or select the desired result record and click the **<Details>** button.

When working with the results, you can perform various actions, such as deleting, reanalyzing, uploading data to the LIS, recalculating the results and printing them. To do this, the user needs to check the box next to the desired result or select several results at the same time by checking the box in the header of the first column and click the button located under the list of results.

Click the **<Search>** button on the **[Result]** screen to go to the **[Result Search]** screen. Then specify the search conditions in the criteria field and click **<OK>**. The list of analysis results displays information about the analysis results from the archive.

The screenshot shows a software interface for viewing analysis results. At the top, there are tabs for 'Current database' and 'Filter', along with search and quick search buttons. A table displays the results, with columns for various parameters. The first row of data is highlighted. Below the table, there is a large empty space for additional information or actions. At the bottom, a toolbar contains buttons for 'Information', 'Delete', 'Repeat', 'Recalculati', 'Network', 'Export', and 'Print'.

ID of dilution samples	Status	Measurement	Concentration	Concentration CV	Label	Repeats	Normal range	Measureme. Module	Load	No. in LIS
101	FT4	41.3 pg/nl 75000	0	0	End	3	2022-02-23 14:02:31	XS_Q2		

Fig. 4.5-1. Results

4.6 Completion of the analysis

Upon completion of all necessary analyses, the user must log out of the user account. Software and Windows operating system before turning off the power to each component.

4.6.1 Preparing to put the analyzer into sleep mode or turn it off

After the analysis is completed, it is necessary to prepare the analyzer for hibernation or shutdown.

1. Remove all samples, quality control solution from the sample area.
2. Remove all reagents from the reagent area.
3. Empty the liquid waste tank.
4. Empty the containers for used cuvettes and tips, and dispose of the blocks of used cuvettes and tips.

4.6.2 Power off

After completing all the analyses, the user can log out of the user account. Software and disable the computer, analyzer and other system components, or it can put the analyzer into sleep mode by selecting the desired option in the software.

1. When the analyzer is in sleep mode, power is supplied only to the reagent area cooling system, and power to the other modules is turned off.
Set them to sleep mode ([Home screen] > <Automatic Sleep Mode Switching >).
2. Turn off power to the analyzer and all modules.

Note.

If the analyzer is suspended for an extended period of time, turn the chemiluminescence immunoassay module and sample processing system switches to the "off" position (on the rear control panel).

4.6.3 Actions after power failure

Check the following in preparation for subsequent assays:

1. Check for unused cuvettes on the horizontal surfaces.
2. If the manual preparation mode is set, a solution of flushing liquid in working concentration must be prepared in advance.
3. Check the flushing fluid tray and starter reagent container for dirt and stains. If present, remove dirt and stains, wipe with a clean soft cloth.
3. When using the liquid waste tank located at a distance, it is necessary to empty the liquid waste into the liquid waste tanks located at a distance in each module;
4. Check the flat part of the analyzer for dirt. If present, remove dirt and stains, wipe with a clean soft cloth.

4.6.4 Removal of the reagent kit

1. Open the door of the reagent area.
2. Holding the handle, remove the reagent kit from the reagent area. If there are no tests left in the reagent kit, dispose of it properly. If reagents are present in the kit, cover the cells of the kit with additional sealing film and place the kit in a special container provided by our company to avoid liquid evaporation.
3. Place the reagent kit in the refrigerator, positioning it appropriately.

4.6.5 Removing and maintaining the sample rack

4.6.5.1 Removing the sample rack

1. Open the door of the sample area.
2. Holding the handle, gently pull the sample rack and remove it completely from the sample area.
3. If there are no samples in the tubes, dispose of the tubes according to the relevant regulatory requirements. If the tubes contain samples that can be used for further analyses, seal them and store them according to the requirements established by the laboratory.

4.6.5.2 Sample rack maintenance

In certain cases, maintenance of the sample rack may be required.

- **Case 1:** after direct contact with the sample.
- **Case 2:** sample tubes cannot be secured in the sample rack.

In case 1, perform the following steps.

1. Transfer all samples from the contaminated rack to a clean sample rack.
2. Remove the sample rack and wipe it dry with a clean, dry towel. To disinfect the rack, spray the contaminated areas with ethanol.

In case 2, perform the following steps.

1. Take a pointed instrument, pry the tip of the instrument from the inner cell retainer of the sample rack as shown in Figure 4.6-1, and gently lift the retainer.
2. Insert the assay tube with the appropriate specifications and check that it is securely fastened.
3. If the tube cannot be secured, follow steps 1 and 2 until the tube is secure.


Fig. 4.6-1. Lifting the retainer

4.6.6 Waste disposal

4.6.6.1 Disposal of used cuvettes and disposable tips

● Disposal of used cuvettes

Used cuvettes are placed in the used cuvette container located on the lower right side of the analyzer housing. The software automatically counts the number of spent cuvettes. You can view the statuses on the [Home Screen] and the main page of the chemiluminescence immunoassay module.

Slide out the cuvette container, remove the waste bag and replace it with a new one, then close the container. The removed waste bag must be sealed tightly and then disposed of properly.

Warning

- 1) After emptying the used cuvette container, reset the number of spent cuvettes to zero.
- 2) Since substances left in the cuvettes can lead to infection, the cuvettes should be placed in a waste bag. Used cuvetts and liquid waste should be disposed of in accordance with local statutory and regulatory requirements.

● Disposal of used disposable tips

Used disposable tips are placed in the tip container located in the middle left of the analyzer housing. The number of spent disposable tips and status can be viewed on the [Home Screen] and the chemiluminescence immunoassay module home page.

Button	Indicator light color	Description
Tip container switch	Green	The tip container is on standby and ready for use.

		When the button is pressed, the sliding compartment will open or close by itself
	Flashing green	The tip container is in the startup state. Pressing the button will stop moving the used tip container immediately, pressing the button again will continue moving the container.
	Orange	The tip container is in the "occupied" state, the button is not available.
	Flashing orange	The tip container error, button not available. In this case, push (or pull) the tip container so that it starts to move inward (or outward).
	Doesn't light up	The analyzer has not been initialized, the tip container can be opened or closed manually.

The procedure of actions%

1. Log in to the user software. Once the analyzer initialization is complete, the number of spent disposable tips can be viewed in the software.
2. When the used tip container button light turns green, click the used tip container button; the light will then turn orange. When the analyzer completes the current operation, the used tip container will begin to slide out of the analyzer and the button indicator will flash green.
3. Wait until the disposable tip container button light stops flashing and is solid green, then empty the disposable tip container.
4. After emptying the container, click the used tip container button; the button indicator light will flash green. Wait until the used tip container is moved inside the analyzer and in the proper position; then the button light will turn solid green. This will indicate that the procedure for emptying the container for recyclable disposable tips has been completed.

Warning

- 1) After emptying the used tip container, reset the number of used disposable tips.
- 2) Since substances left on disposable tips can cause infection, the tips should be placed in a waste bag. Used disposable tips and liquid waste should be disposed of in accordance with local legal and regulatory requirements.

The warning. The container for used cuvettes has the following dimensions: 466 x 219 x 204 mm (length x width x height); the dimensions of the container for used tips are: 374 x 118 x 200 mm (length x width x height). The user can purchase a suitable waste bag that complies with local regulatory requirements.

- Disposal of used cuvette and tip assemblies

Used cuvette blocks are placed in the sliding cuvette storage compartment, and used tip blocks are placed in the sliding tip storage compartment. Refer to section 4.2.4.1 "Loading Cuvettes and Disposable Tips" for retrieval information. Used cuvette and tip assemblies do not pose a risk of biological contamination and should be disposed of following standard waste disposal procedures.

4.6.6.2 Liquid waste disposal

The outlets for high and low concentration liquid waste are located at the bottom of the rear panel of the analyzer and are labeled accordingly.

The liquid waste can be stored in a dedicated liquid waste tank or disposed of in accordance with local legal and regulatory requirements. When using the dedicated liquid waste tank, you can view the statuses on the home screen and home page of the chemiluminescence immunoassay module. If it is indicated that the waste tank is full, waste should be disposed of in a timely manner.

Warning

- 1) Always wear gloves when disposing of waste.
 - 2) The liquid waste should be disposed of in accordance with local legal and regulatory requirements.
 - 3) Clean the liquid waste tank regularly according to the maintenance documents.
-

Chapter 5 Software functions

5.1 Description of the software interface

After logging into Windows, double-click the user software shortcut on your desktop. After starting the software, a login dialog box will appear on the screen. Enter your user name and password, and then click <Log in> to go to the home screen of the user software. The description of the software interface is given below.


Fig. 0-1. User software interface

5.1.1 Main menu

The main menu contains eight buttons.

Fig. 0-2. Main menu of user software

The buttons have the following functions.

Button	Function
Sample	Registration and viewing samples, performing other activities.
Result	Check or delete sample results, perform repeat measurements, and other operations.
Reagent	Check reagent information, record reagent calibration, check calibration results, or perform other operations.
QC	Registration and viewing of QC results, calculating statistics on QC results, viewing a QC graph and performing other operations.
Utility	Registration of data for system verification, viewing the results of this verification, as well as performing other similar operations.
Settings	Configure analysis settings, profile, network, language, etc.
Home Screen	This screen displays the current status of analyzer components, analysis time countdown, etc.
Tool icon	Shows the current status of the analyzer.

5.1.2 Interfaces

Displays details of each main menu function and is named according to those functions. The figure below shows the **[Sample]** screen.


Fig. 0-3. Screen [Sample]

5.1.3 Submenu

Select **<Settings>** from the main menu to access the submenus of this function. The name of the submenu is the same as the name of the buttons in the main menu. The figure below shows the **[Settings]** menu.


Fig. 0-4. Menu [Settings]

5.1.4 Toolbar

The toolbar contains several shortcut buttons and an area with statuses.

1. The shortcut buttons are used to quickly start a specific function.

- : **<Warning message>** button. When an alarm is generated by the analyzer, this button indicates the type of error that occurred. By pressing this button, the user can view the details of the alarm.
- : **<Start>** button. When this button is pressed, upon completion of sample registration, calibration, and control material, the software will send a command to the analyzer to start and run tests.
- : button **<Suspend transport >**. Click this button if you want to pause the sample racks data output in the sample area during the analysis run. The analyzer will perform this action and the button will change to **<Resume Transport.>**. Click the button **<Resume Transport.>** to perform this operation.
- : button **<Stop analyzer>**. If you need to stop the analysis, click this button, select the module you want to interrupt in the dialog box that appears and click **<OK>**. The user software will send a stop command to the analyzer and the analysis in the selected module will immediately stop and its status will change to "Stopped".

Note.

1. Use of this function is not recommended unless special circumstances warrant, such as analyzer malfunction. Once the assay is stopped, all incomplete assay procedures will become invalid.

- : button **<Exit>**. Press this button if you want to display the [Off/Exit Account] screen where you can lock the screen or exit the software. When you press **<Lock Screen>**, the lock screen is displayed, where you can change the user account or lock the software. By entering the user name and password, you can log back into the software. In the locked state, data exchange between the user software and the analyzer will proceed as normal without stopping the analysis procedures. The screen will display the current system time, number of samples and assays, and error messages. Click on the light indicator icon in the lower right corner of the screen if you want to disable the alarm sound speaker.

Click the **<Exit>** button to exit the control software.


Fig. 0-5. Lock screen

- : button **<Review Details>**. Touch this button to display the [Review Details] screen. In the dialog box that opens, you can view all assays registered for each sample, the number of registrations for each assay, and the number of remaining tests for the reagent.


Fig. 0-6. Review Details

-  - LIS status indicator (LIS connection status). The indicator lights up white when the connection to the LIS server is successful and red when the connection fails.

5.1.5 Dialog boxes

Dialog boxes open when you press a specific button on the screen. You can open one or more subwindows in the dialog box. The dialog boxes are divided into function windows and information confirmation boxes.

5.1.6 Buttons

Menus, screens, and dialog boxes contain one or more buttons that are used to launch functions. The three most common buttons are **<OK>**, **<Cancel>** and **<Back>**. The **<OK>** button in the dialog box is used to confirm the operation; the **<Cancel>** button in the dialog box is used to cancel the operation; the **<Back>** button is used to close the dialog box without performing the operation specified in the message of this window.

5.2 Function description [Sample]

Click the **<Sample>** button in the main menu to go to the **[Sample]** screen. The functions of the **[Sample]** window are used to register tests, check the test list and other operations.

5.2.1 Sample registration






Fig. 5.2-1. Sample

5.2.1.1 Decoding of symbols

Area	Symbols	Description
Sample area	1	The sample rack is not loaded
	2 A002	The rack with unscanned samples is loaded
	6 NA006	The standard sample rack is loaded
	7 CA007	The tripod for the QC is loaded
	2 EA002	The STAT sample rack is loaded
	4 LA001	Tripod loaded to test the light signal
	5 NA005	The patient sample rack loaded
	8 NA008	Selected sample rack
	10 NA010 ✓	Massively selected tripods with samples
	1 ⚠ NA001	You cannot remove a tripod with samples in the testing state
	12 ERROR	The number of the sample tripod cannot be identified or duplicated

Analysis selection area		The set of reagents for this analysis is installed in the reagent area, the remaining number of tests is greater than zero, the analysis is assigned to the analyzer.
		The set of reagents for this analysis is installed in the reagent area, the remaining number of tests is zero, the analysis is assigned to the analyzer.
		The reagent kit for this analysis is not installed in the reagent area, and the analysis is not assigned to the analyzer.
		Analysis profile

5.2.1.2 Registration methods

1. Registration of a single sample

Load the samples for analysis into the sample stand. If there is a barcode on the test tube, turn the test tube with the barcode to the open side of the tripod. Load the tripod into the sample area properly. Click the **<Scan>** button (see the section 5.2.2 "Scanning"). Select the sample rack whose information is to be changed and enter the sample IDs in the tube positions one at a time (not required if a barcode is used). Select the sample ID and assay from the assay selection area on the right side of the **[Sample]** screen or select an assay profile (if you want to quickly select multiple assays) to complete the sample registration.

Registration of STAT samples is performed in the same way as registration of standard samples, however, STAT tripods are given priority when performing analysis.

If the sample needs to be diluted, click **<Dilution>**, select analysis in the **[Dilution Coefficient] dialog box.**, specify the coefficient and click **<OK>** to register the dilution coefficient.

When using the LIS function, you only need to set the settings correctly and connect the LIS server. After entering the sample data, the analysis data for the sample will be automatically received by the analyzer from the LIS server. Another way: Click **<Use for LIS>** (see 5.2.3 "Use for Analyses") to manually select a sample and get the analysis data for the sample from the LIS server.

2. Registration of the series of samples

Select one or more sample racks in the sample area and specify the ID of the first and last sample. Then press the Enter key to register a series of the ID samples on multiple racks. Then select the analysis or use the profile function to select the analyses to register the samples.

For information about recording diluted samples or STAT samples, as well as information about using the LIS function to obtain analysis data, see in section 1, "Registration of a single sample".

Note:

The priority of tripods for transporting samples is as follows: tripod for checking the light signal > tripod with STAT samples > tripod for patient samples > tripod for

quality control > tripod with standard samples.

5.2.2 Scanning

After loading the tripod with samples into the sample area, the analyzer first scans the tripod to identify it. Click **<Scan>** to go to the **[Sample Tripod Scan]** screen, where you can specify which tripods to scan. Click **<Scan>** and the analyzer will start scanning the selected tripods in turn. The sample processing status will change to "Scanning". The samples cannot be recorded during the scan.

On the **[Sample Tripod Scan]** screen, you can click the button in the "Outpatient" column of the standard sample tripod and mark it as the patient's sample tripod. The value in the column will change to "Yes".


Fig. 5.2-2. Scanning a tripod for samples

5.2.3 Application for analyses

Tap **<Use for LIS>** on the **[Sample]** screen to go to the **[Use for LIS]** screen. On the **[Usage for LIS]** screen, you can view the status of the sample in the LIS system. To get the analysis data from the LIS server and see the progress of the analysis in LIS, select the sample and click **<Use for LIS>**. To stop using the analysis in the LIS, click **<Cancel>**.


Fig. 5.2-3. Use for LIS

5.2.4 Worklist

On the [Sample] screen, press <Worklist> to go to the [Work list] screen. The functions include displaying a list of registered sample analyses, canceling analyses for which samples have not been added, changing the number of analysis repeats, deleting already registered sample analyses, as well as exporting information from the sample analysis list in Excel format.


Fig. 5.2-4. Worklist

1. Cancel analyses without uploading samples

If, after completing the sample registration and starting the analyses, it is necessary to cancel the test, the user can click the <Cancel Test> button to go to the screen [Cancel Test]. The list displayed on the screen will contain the analyses that can be canceled at the moment (the samples for analysis are still in the sample area). Select the analyses that you want to cancel (filtering by tripod number and analysis is available), and then click the <Cancel Test> button to cancel the selected analyses. Since some tests may not meet the cancellation requirements during the procedure, the list for tests will indicate which tests were actually canceled.


Fig. 5.2-5. Cancel analysis

2. Change the number of repeats

You can change the number of repetitions of one analysis using the "+" and "-" buttons in the column indicating the number of repeats. If you want to change the number of repeats for a series of registered

tests, select the desired tests, click the <Repeat Number> button to access the [Change Repeat Number] screen, and then change the number of repeats with the + or - buttons or manually enter the desired number in the "Change Value" column.


Fig.5.2-6 Change the number of repeats

3. Deleting the analysis

Select the sample analysis that you want to delete, click the <Delete> button and confirm the action, only registered analyses can be deleted.

4. Exporting analysis data

Select the analysis of the sample whose information you want to export, and click the <Export> button to export the information in Excel format.

5.3 Description of the screen functions [Result]

Click the <Result> button in the main menu to go to the [Result] screen. The [Result] screen is used to view, search and delete results, view detailed information about the results, repeat measurements and upload data to the LIS, as well as to recalculate, export and print the results. This section describes the functions of the [Result] menu sequentially.


Fig. 5.3-1. Results

5.3.1 Search for results

On the **[Result]** screen, click the **<Search>** button and set the search criteria on the screen **[Search for a result]** so that the required results are automatically filtered. To display the results of all analyses for the current day, click the "Today" button. The "Results" field in the upper right corner of the screen displays the total number of entries in the current list.

Enter the sample or assay ID in the quick search box or click  to the right of the edit box to automatically scan the assay results. The results of the analyses corresponding to the selected criteria will be displayed as a list.

After transferring the backup (see Section 5.7.2.8 Backup copy/View), click **<Current>**. The button will change to **<Log>** button, and the imported history data (backup copy) matching the search criteria will be displayed in the search results list.

5.3.2 Deleting the results

From the list of test results on the **[Result]** screen, select the result you want to delete and click **<Delete>**.

Note.

1. Only the results of tests that have been completed (including those with an error) can be deleted.
- 2 Please note that deleted test results cannot be recovered.

5.3.3 Information about the result

From the list of test results on the **[Result]** screen, select a result and double-click on it, or click **<Information>** to go to the **[Result Details]** screen and view detailed information.

№	RLU	Concentration	Status	Label	Final cycle	End time	Label
1	75000	41.3	end	>	136	2022-02-23 14:02:31	> More than the upper limit of the norm

Fig. 5.3-2. Information about the result

Press <Convert to QC> to go to the [Convert Sample Data to QC Data] screen. Select the batch number for the QC in the Batch No. field and click <Save>. The results of the analysis will be transferred by the software to the results for the current day for the corresponding batch of control materials.


Fig. 5.3-3. Converting sample data to QC data

Click the <Export> button to export the result details in Excel format.

5.3.4 Repeat measurement

From the list of test results on the [Result] screen, select the test results and click the <Repeat Measurement> button. After successful registration of the repeat measurement procedure, the following message will appear on the screen: "Registration successful". In the list of assay results to which sample remeasurement results have been added, check the remeasurement results marked with an asterisk.

Note.

1. It is not possible to remeasure the registered samples.
2. The reanalysis sample rack must be in the sample area, it cannot be removed.
3. The remeasurement parameters must match those of the original analysis. If you change analysis parameters, you must register it again.

5.3.5 Sending results

After enabling the LIS function and successfully connecting to the LIS server, mark the test results to be sent in the [Result] window and click the <Network Connection> button to upload the test results to the LIS server.

5.3.6 Recalculation of results

If the analyte concentration is to be recalculated based on a new calibration curve, select the analyte to be recalculated from the list of analyte results on the [Result] screen and click <Recalculate> to start the recalculation. After recalculation of the result, the value obtained will be displayed for the corresponding assay in the list of assay results.

Note.

After using the result re-calculation function, the original concentration value cannot be restored. Be careful!

5.3.7 Exporting results

In the list of analysis results on the **[Result]** screen, select the analysis results that you want to export and click the **<Export>** button to export the results in Excel format.

5.3.8 Printing the results

In the list of test results on the **[Results]** screen, select the test results that you want to print and click the **<Print>** button to go to the **[Preview]**. Here you can view the results that you want to print. Check if the information is correct and click on the print icon in the upper-left corner of the print dialog box (an external printer must be connected).


Fig. 5.3-4. Preview

5.3.9 Description of result labels

On the **[Results]** screen, the results list will display the tag information of each result to help users evaluate the reliability of the results. If an anomaly is detected as a result, repeat the experiment after the anomaly has been eliminated. Information about the label of the immune module of this analyzer is given below.

Table. 5.3-1. Description of result labels

Label	Description
S	Sample aspiration error
L	The sample probe does not detect liquid
K	Not enough reagent

Label	Description
G	The left reagent dispenser does not detect the liquid
I	The right reagent dispenser does not detect the liquid
D	The probe of the samples found a clot
T	Device malfunction
N	Sample probe error
M	Reagent probe error
*	Analyzer failure
E	The reagent has expired
C	Calibration is invalid
R	Recalculation of results
<Q	Less than the lower limit of the QC range
>Q	More than the upper limit of the QC range
<<	Below the lower limit of the measurement range
>>	Above the upper limit of the measurement range
<	Below the normal range
>	Above the normal range

5.4 Function description of the [Reagent] screen

Click the <Reagent> button in the main menu to go to the [Reagent] screen. The [Reagent] screen functions are used to view reagent and calibration information, log calibration data, and view calibration results and reagent consumption information. This section consistently describes all the functions of the [Reagent] screen.

5.4.1 Information about the reagent

On the [Reagent] screen, the user can view information about the right and left areas of the reagents, such as information about the installed reagent kit, status, name, and mixing time of magnetic microparticles.



Fig. 5.4-1. Screen [Reagent]

5.4.1.1 Reagent color marking

Symbols	Description
	The number "1" indicates that the reagent is set to position 1 in the reagent area, and the background color indicates the calibration of this reagent.. Each calibration state has its own color marking:  — valid calibration data  —there is no valid calibration data.  —new calibration data.  — calibration is in progress  —the validity period of the calibration expired

	"FT4" indicates the name of the current reagent (if a reagent loading error occurs, the message "ERROR" will be displayed here); "95" indicates the number of remaining tests for the reagent; the icon  indicates that the reagent is currently in use and the reagent kit cannot be extracted; the background color indicates that the expiration date the reagent has not expired. Each reagent state has its own color marking:  — the reagent has expired  — the reagent is not installed in the specified position  — normal condition, reagent not selected  — the position is selected
	 means the current reagent is disabled and this reagent kit cannot be used for analysis.
	"00:00" — countdown of the mixing time of magnetic microparticles of the current reagent
	A check mark indicates that the current reagent is selected

5.4.1.2 Description of the time state

By default, the mixing time of magnetic microparticles is 30 minutes. If the user clicks the **<Start>** button on the action bar when there are less than 00 minutes left before the end of the mixing process, a warning will appear on the screen that the magnetic microparticles are not completely mixed. The analysis can be started only after the user clicks the **<OK>** button.

Warning

If the magnetic microparticles are not completely mixed within 30 minutes, this may affect the results of the analysis.

5.4.1.3 The procedure for using reagents

More than one reagent kit can be loaded for each reagent analysis. In this case, the reagents are used in the following order:

- 1) Reagents with the lowest residual volume should be used first.

2) Reagents are used from left to right.

5.4.2 Curve query

On the **[Reagent]** screen, click **<Curve Query>** to go to the **[Curve Information Query]** screen, and select Module and Analysis. All calibration curves of the corresponding reagent analysis will be shown in the list in the lower left area. In addition, you can select a calibration curve from the list on the right side of the screen to view detailed information about the curve. To print information about the selected calibration curve from the list, click **<Print>**.


Fig. 5.4-2. Curve query

5.4.3 Calibration registration

On the **[Reagent]** screen, select one or more reagents to calibrate and press **<Reg. calibr.>**. In the dialog box that appears, confirm the selection of the analysis for calibration and click **<OK>** to register the calibration analysis. See the section 4.3.1 "Registration of calibration analysis".

5.4.4 Calibration results

On the **[Reagent]** screen, click the **<Res. calibr.>** to go to the **[Res. calibr.]**. The calibration results for all reagents will be displayed in the list at the top of the screen. When you select a calibration result, the information corresponding to this result will be displayed in a list at the bottom of the screen. All the same operations are available for the results of the system check as for the results of the analysis. For more information, see in the section 5.3 "Description of the functions of the [Result] screen".

Fig. 5.4-3. Calibration results

5.4.5 Information about the reagent

In the [Reagent] window, double-click on the reagent or select it and click the **<Calibr.> button. the curve.>** to open the [Reagent Information] screen. The dialog box will display basic information, curve and calibration information, and the corresponding functions for the reagent will be available.

Fig. 5.4-4. Information about the reagent

The following is a description of the parameters from the calibration information area.

Content	Description
Calibrator 1, calibrator 2	The results obtained for two calibrators
RLU	The measured RLU value of the calibrators (up to three measurements for each calibrator)
Average RLU	The average value of the three dimensions of the RLU
RLU – CV	Coefficient of variation for three dimensions of the RLU
RLU – deviation.[%]	The difference (in percent) between the average value of the three RLU measurements and the factory standard value
Expected concentration	The set calibration concentration value

Calculated concentration	Click the <Calculate Curve> button to calculate the concentration for the last or new calibration curve
Concentration deviation (%)	The difference (in percentage) between the expected concentration and the actual calibration value
The range of permissible concentration deviations	The permissible range of deviations of the measured calibration concentration values
Permissible deviations of the RLU	The permissible range of deviations of the RLU values during calibration.

5.4.5.1 Operations performed with calibration results

When the calibration analysis is complete, open the **[Reagent Information]** screen to view calibration and operating curve information. If the calibration result exceeds the acceptable limits, the corresponding calibration parameters will be highlighted in red. The user can evaluate the results and perform the following actions.

1. When you receive a valid calibration result, click the **Confirm** button, and then the **<Confirm>** button. The calibration result will be accepted by the system.
2. If you receive an invalid calibration result, click **Reject**, and then click **<Confirm>**. The system will reject and delete the calibration results. After that, the calibration will need to be performed again.
3. If you have selected the **Recalculate function for related samples** or **Recalculate samples with the same kit number** and clicked **<Confirm>**, the latest calibration curve will be used to recalculate the concentration for all the results of the current analysis (applicable only to the results of analyses of samples and control materials for the current day).
4. To print the current calibration information, click the **<Print>** button.
5. You can change the concentration tolerance and the RLU tolerance.

Note.

1. It is necessary to adjust the permissible concentration deviation and the permissible RLU deviation carefully, since this may affect the reliability of the analysis result.

5.4.5.2 Calibration log

On the **[Reagent Information]** screen, click the **<Calibration Log>** button to go to the **[Calibration Log]** screen.

Calibration log					
Module :	XS.02		Reagent	F14	
batch no.:	2021010521010		kit no.:	2501	
			Reagent area	Left reagent	
			ID reference curve:	2021012113030	
Date calibr.	Calibrator1 - Average	Calibrator2 - Average	Calibrator1 - CV(%)	Calibrator2 - CV(%)	Calibrator1 - deviation
2022-02-23 16:04:14	581251	169293	1.52	0.17	-5.6

Fig. 5.4-5. Calibration log

Select an entry in the calibration log and click the **<Restore Calibration Data>** button to replace the current reagent calibration data with the selected data. Click the **<Delete>** button to delete an entry from the calibration log.

In the calibration log, double-click on an entry or select it and click the **<Details>** button to go to the **[Calibration Log Details]** screen to view calibration information.

To print the current calibration log, click the **<Print>** button.


Fig. 5.4-6. Calibration Log information

5.4.6 Reagent Information List

On the **[Reagent]** screen, click the **<Reagent Information List>** button about the reagent> to open the **[Reagent Information List]** option. The list of all reagents that have been loaded into the analyzer will be displayed.

To delete reagent kit entries, select one or more entries and click the **<Delete>** button, after which the selected entries will be deleted. To export records, click **<Export>**, after which the selected records will be exported in Excel format. To print reagent records, press the **<Print>** button, after which the selected records will be printed.

Nr	Analysis	Module	batch no.	kit no.	Remaining tests	Open (days)	Suitability
1	F13	K1_02	202105170222	6	100	1	2022-01-09

Fig. 5.4-7. Reagent Information List

5.4.7 Reagent consumption statistics

On the [Reagent] screen, click the **<Reagent Consumption>** button to open the [Reagent Consumption Statistics] screen. Statistics on the consumption of reagents used to perform analyses can be viewed for a certain period of time.

Select the start and end date and click the **<Search>** button, after which the screen will display up-to-date data on the consumption of reagents for analysis for the specified period of time.

Select one or more information records and click the **<Export>** button to export the selected information in Excel format. To print the selected reagent consumption statistics, click the **<Print>** button.

Nr	Type reagent	Name reagent	Sample (number of times)	QC (number of times)	Calibr. (number of times)	Total (number of times)
1	Immun.	Ice				

Fig. 5.4-8. Reagent consumption statistics

5.4.8 Enabling and disabling the use of reagents

Select the unconnected reagent kit and press **<Off/on>**. The reagent status will change to "Disabled" and the analyzer will not use it for analysis. Select the reagent kit, press **<Off/on>** again, and the kit status will change to "Enabled".

5.5 Description of the screen functions [QC]

Click the <QC> button in the main menu to go to the [QC] screen. The functions of the [QC] screen are used to register QC data, work with the results QC and viewing the QC charts of the batch and QC data for the month. This section describes all the functions of the [QC] screen sequentially.

5.5.1 QC registration

On the [QC] screen, click on the <Reg. of data for QC> tab to go to the [Reg. of data for QC] screen. It is designed to display information about materials for QC and QC analyses, as well as adding, changing and deleting analyses and QC materials.



Fig. 5.5-1. Registration of data for QC

5.5.1.1 Imp./exp.

Press <Imp./exp.> to go to the [Import/Export of QC Materials] screen. After selecting the path type, you can export the material settings for QC from the software to the specified path or import the material settings for QC from the specified path to the software.

5.5.1.2 Adding materials for QC

Click <Add> on the <Reg. of data for QC> screen to go to the [QC Material Settings] screen. Enter the name, lot number, tripod position, level, manufacturer, expiration date, barcode, module and ID of the sample material for QC. Click <Save>.

Material settings for QC

Name	
batch no.	
Tripod	CA001-1
Level	1
Manufacturer	
Expiration date	2029-02-28 12
Barcode	
Module	300
Sample ID	

Save Cancel


Fig. 5.5-2. Material settings for QC

5.5.1.3 Changing the material data for QC

On the [Reg. of data for QC] screen, select the material entry for the QC and press <Edit> to go to the [Material Settings for QC] screen. On it, you can change the name, lot number, tripod position, level, manufacturer, expiration date, barcode, module and ID of an existing QC material sample.

5.5.1.4 Editing information about the QC analysis

On the [Reg. of data for QC] screen, select the recording of the material for the QC and press <QC Analysis> to go to the [Setup] screen. QC analysis]. On it, you can add analyses of the material for QC and change information about existing QC analyses.

QC analysis setup

☐ Analysis	Number of	Target value	Expected	Expected	User

+ Add Save Delete

Recalculation About QC

Name QC: QC3

batch no. for QC: 1002

Analysis: →

Number of repeats: 1

Target value:

Expected value:

Waiting CV:

Expected range: -

User: -

Return

Fig. 5.5-3. QC Analysis Settings

5.5.1.5 QC registration

On the [Reg. of data for QC] screen, select one or more material entries for the QC, and then select the desired module, tripod number, and position. Click <Registration> to register the material.

5.5.2 QC results

Click on the <QC Result> tab on the [QC] screen to go to the [QC Result] screen. This screen is used to search and delete QC results, view detailed information about QC results, send QC results to the LIS, as well as to recalculate, export and print QC results. All the same operations are available for the results of the system check as for the results of the analysis. For more information, see in the section 5.3 "Description of the functions of the [Result] screen". For a description of the result labels, see Table 5.3 1.


Fig. 5.5-4. QC results

5.5.3 Analysis of control materials

To view and print QC graphs, click on the <QC Analysis> tab on the [QC] screen and go to the [QC Analysis] screen.


Fig. 5.5-5. QC analysis

5.5.3.1 Viewing QC graphs

- QC graphs for the party

Remove the selection from the **Monthly QC** on the left side of the **[QC Analysis]** screen and specify the search criteria. The search results will be displayed in a list in the lower right part of the screen. Select an entry in the list, and then the corresponding curve for the selected batch will be displayed in the lower area of the QC graph window.

- QC graphs for the month

If you need to check the QC results for one month, when specifying the search criteria, select the **"Monthly QC"** checkbox, this will allow you to evaluate the QC results for one month and build an appropriate QC graph. For further steps, see **"QC graphs for the batch"** above.

5.5.3.2 Operations performed with the results of the QC

In the QC analysis tab, in the lower-left corner of the QC project list, select one control element, then mark the required batch number of the control to view the corresponding quality control results, click the **<QC Data>** button. The **[Ignore QC data]** screen opens with a list of detailed QC results corresponding to the batch number of control materials. In the list, select the QC results that you want to ignore. They will be ignored, and the graph The QC and the list of conclusions on the QC have been compiled anew. In the monthly QC mode, click the **[Ignore QC data]** button in the **<Ignore QC data>** window to ignore all points except the last one for the current day.

Select the QC entry in the list of QC results located in the lower left part of the screen, after which a list with QC conclusions will be displayed in the upper right part of the screen. The QC conclusions are evaluated in accordance with the QC user rules. For more information, see the section 5.7.1.1 "Analysis Parameters".

5.6 Description of the screen functions [Utility]

To go to the **[Utility]** screen, click the **<Utility>** button in the main menu. On-screen functions **[Utility]** is used to register BGW and LC tests, as well as to view the results of these tests in the immunochemiluminescence analysis module. This section discusses the functions presented on the screen **[Utility]**.

5.6.1 System checks

5.6.1.1 Registration of system checks

On the screen **[Utility]** click **<System checks>** under the **"Immunn."** panel on the left side to go to the **[System checks]** screen. Then click **<Add>** to open the **[Registration of data for system verification]**

screen, as shown in the figure below.

On the **[Registration of data for system verification]** screen, enter the numbers of BGW checks, LC checks of the left and right reagent dispenser, and LC checks of the sampler, and select the positions of the LC sampler tripod. Then click **<Registration>** and **<OK>** to complete the registration of the system verification data and start the process. For more information about the system check, see in the section 4.2.4.2 "Confirmation of the results of the system check".

Fig. 5.6-1. Registration of data for system verification

Note.

1. Before performing the LC test, place the light signal test solution in the appropriate positions in the sample area, as well as on the right and left reagent areas. For more information, see the section 5.6.2 "Loading liquid to check the light signal".
2. BGW and LC checks must be carried out after the replacement of the starting reagent. For more information, see the section 4.2.3.2 "Adding flushing liquid and starting reagent".

5.6.1.2 System check results

On the screen **[System checks]** you can search, delete, view information, export, and print system check results. All the same operations are available for the results of the system check as for the results of the analysis. For more information, see in the section 5.3 "Description of the functions of the [Result] screen".


Fig. 5.6-2. System check

**Warning**

The analysis is allowed only upon receipt of the results and coefficients of variation of the system checks that meet the requirements. It is not recommended to perform sample analysis if the test results are outside the acceptable range (deviating values are highlighted in red in the system test results list). Otherwise, it may adversely affect the validity of the assay results.

5.6.2 Loading liquid to check the light signal

The Light Signal Check (LC) is used to check the dosing accuracy of the sampler holder and the reagent dispenser holder of the analyzer. When performing the LC check, observe the following requirements:

- The liquid for testing the light signal must be placed in the first tube of the LC sample stand and the tripod must be installed in the sample area.
- The LC test liquid for the left reagent dispenser is placed in the 4th hole of the 21st cell of the left area, and for the right dispenser - in the 4th hole of the 42nd cell of the right area.

Strictly follow the instructions for preparing and storing the liquid to check the light signal.

5.7 Description of the functions of the [Settings] screen

Click the <Settings> button in the main menu to go to the [Settings] screen. The functions of the [Settings] screen are used to configure various parameters and set system values, as well as to perform analyzer maintenance and other actions. This section describes all the functions of the [Settings] screen sequentially.


Fig. 5.7-1. Settings

5.7.1 Setting Parameters

The first column of the [Settings] screen is [Analysis parameters]. It allows you to change analysis parameters, profile definitions, key positions, result format, result definitions, purpose of analysis, preset volumes of reagents and warnings about reagents. The following provides basic information about the various parameter settings.

5.7.1.1 Analysis Parameters

Click the <Analysis Parameters> button to go to the [Analysis Parameters] screen shown in the figure below. Configure the appropriate parameters for each analysis. On the analysis parameters screen, you can set the parameters of the chemiluminescence immunoassay module. The parameters of the analyses are given below.

Select the Analysis information from the list on the left side of the <Analysis Parameters> screen, and on the right side there are the following tabs: basic parameters, reference curve, QC rules, dilution coefficient, automatic dilution and reflex testing.

1. Basic parameters

Select the analysis information from the list on the left side of the [Analysis Parameters] screen, click on the <Basic parameters> tab.. The [Basic parameters] screen will appear] with information about the analysis.

1) Changing the analysis information

Select the analysis data from the list and change the required information on the [Basic parameters] screen. Then click <Save>.

2) Deleting the analysis

Select the analysis data that you want to delete from the list and click <Delete>.

Fig. 5.7-2. Analysis Parameters

Note.

You can delete the analysis data from the list, however, deleting the analysis information may affect the original information about it. Delete information with care.

2. Reference curve

Click on the **<Reference Curve>** tab and select the analysis data from the list on the left. The corresponding ID value will appear in the "Reference Curve ID" area of the Reference Curve screen, as shown in the figure below. The reference curve is the basis of calibration. For each analysis, there is a reference curve that is assigned a specific identification number and batch number.

Fig. 5.7-3. Reference curve

1) Viewing the reference curve

The ID of the reference curve of the corresponding analysis is shown in the area **[ID of the reference curve]**. Click on the ID of the reference curve from the list to display detailed information about it.

2) Selection of the reference curve

Click on the reference curve in the **[ID of the reference curve]**, change the data, click the

<**Curve Selection**> button to re-fit the reference curve and display the new curve on the screen. Then click the <**Save**> button to complete the curve selection and save the data.

- 3) Removing the reference curve
Select the reference curve and click the <**Delete**> button. Click the <**OK**> button to delete the reference curve.
- 4) Printing a reference curve
Select the reference curve and click the <**Print**> button to print the reference curve and the corresponding analysis data.

3. Quality control rules

Click on the <**QC Rules**> tab to go to the [QC Rules] screen shown in the figure below. The software supports the Westgard rules, allowing the user to set the necessary QC rules for analysis.

Select the analysis data from the list, select the QC rule in the [Westgard Rules] section and click <**Save**> to complete the rules setup.


Fig. 5.7-4. Quality control rules

4. Dilution coefficient

Click on the <**Dilution coefficient**> tab. to go to the [Dilution coefficient] screen, as shown in the figure below. For analyses that use a sample diluent, up to 10 dilution coefficients can be set. The dilution coefficient can be adjusted by the user.


Fig. 5.7-5. Dilution coefficient

- 1) Set the dilution coefficient:
Select the analysis from the list, specify the dilution coefficient in the **[Dilution coefficient selection area] window** and click the button without dilution coefficient in the **[Dilution coefficient.]** to complete the dilution coefficient setting for the selected analysis.
- 2) Changing information about the dilution coefficient
Select an analysis from the list of analyses, then select the dilution coefficient in the **[Dilution coefficient.]** and click the **<Edit>** button to go to the **[Dilution Coefficient Settings] screen**.


Fig. 5.7-6. Changing dilution settings

- 3) Removing the dilution coefficient
Select the dilution coefficient in the **[Dilution coefficient.]** and click the **<Delete>** button to delete the dilution coefficient.

5. Automatic dilution

Click on the **<Auto-dilution > tab.** to go to the **[Auto-dilution.]**, shown in the figure below. When using the automatic dilution function, a sample whose result exceeds the automatic dilution threshold is automatically registered in the dilution sample area, after which the analyzer dilutes the liquid and completes the analysis.


Fig. 5.7-7. Automatic dilution

1) Adding an automatic dilution procedure

Select the analysis from the list and select the coefficient in the [Coefficient dilution] section. If dilution is not set for analysis, this area will be empty. For more information about the definition of dilution methods, see the section "Dilution coefficient" above. Click **<Add>** and the new dilution coefficient will be shown at the top of the list.

2) Changing the automatic dilution threshold

Select an analysis from the list of analyses, select the automatic dilution coefficient from the list and enter a critical value, and then click the **<Save>** button to change the threshold value for automatic dilution.

3) Removing the automatic dilution procedure

Select an analysis from the list of analyses, select the automatic dilution coefficient that you want to delete from the list, and click the **<Delete>** button.

6. Reflex test

Click on the **<Reflex>** tab to go to the **[Reflex]** screen shown in the picture below. Reflex testing is an analysis that is triggered if the results of the test being performed are in the range that meets the criteria for performing reflex testing.


Fig. 5.7-8. Reflex testing

1) Adding Reflex Testing

Let's say it is necessary for reflex testing B to be performed if the value of the analysis result A is in the range of 50-150. To do this, select Analysis A from the list, and in the [Select joint analysis to add] area, select reflex testing B. In the reflex testing settings area B, specify the upper and lower limits (for example, the upper limit is 150 and the lower limit is 50), click the **<Add>** button, and then add the relevant information to the list of information about reflex testing.

2) Changing the parameters of reflex testing

Select an analysis from the list of analyses, in the list of information about reflex testing, select the reflex testing whose parameters you want to change, change the upper and lower limits in the appropriate fields and click the **<Save>** button.

3) Removing Reflex Testing

Select an analysis from the list of analyses, in the list of information about reflex testing, select the reflex testing that you want to delete, and click the **<Delete>** button.

5.7.1.2 Profile setup

To go to the **[Profile Setup]** screen (see the picture below), click the **<Profile Setup>** button. The user can create a profile for combining analyses and set key positions for it (see the section 5.7.1.3 "Setting key parameters"). The profile will be displayed in the profile selection area on the **[Sample]** screen. You can assign multiple analyses to a sample with one click, which simplifies the registration process.


Fig. 5.7-9. Profile setup

1. Adding the profile

Select one or more analyses in the selection area in the center of the **[Profile Setup]** screen. In the **[Information] about profile** field Enter the profile name and click **<Add>** to add a profile.

2. Changing the profile

Select a profile from the **[Profile List]** list, change the details in the profile details area on the left, and click **<Save>**.

3. Deleting the profile

Select a profile from the **[Profile List]** list and click the **<Delete>** button.

5.7.1.3 Setting up key positions

Press **<Key Position>** to go to the **[Key Position]** screen, as shown in the picture below. Key positions are used to set positions for displaying analysis data in the analysis selection area on the **[Sample]** screen.


Fig. 5.7-10. Key positions

1. Adding/changing the page

By default, three pages of key position settings are displayed. To add a page, click **<Add Page>**, and to delete it, click **<Delete page>**. You can drag the page name to change the order of the pages, or double-click on its name to rename it.

2. Setting up key analysis positions

By default, the MGL page shows all analyses for tripods inserted into the analyzer. On this page, you can add or clear analyses and profiles. On other pages, you can first select the position on the page where the analysis or profile should be displayed, select the analysis or profile in the corresponding list. To save an analysis or profile position, click **<Save>**.

3. Resetting key positions

Select a position on the page where the analysis or profile is set, and click **<Leave blank>** to reset the key position settings for analysis or profile.

4. Canceling settings

If you have added or deleted pages or set key positions for analyses or profiles, but want to restore the original settings, click **<Cancel>**.

5.7.1.4 Format of the results

To go to the **[Format of the results] screen** (see the picture below), press the button of the same name. The "Format of the results" function is used to determine the number of decimal places in a numeric range.

Specify the upper limit of the numeric range and the corresponding number of decimal places. Specify the upper limit of the numeric range, and the value of the lower limit of the next range will change automatically. Click the **<Save>** button to complete the configuration of the result format settings.

Range	Number of signs
0	3
1	2
10	1
100	0

Save Return

Fig. 5.7-11. Format of the results

5.7.1.5 Setting up results

Click the <Setting up results> button to go to the [Setting up results] screen shown in the figure below. The user can configure the display parameters of the results so that the necessary information is displayed when receiving data.


Fig. 5.7-12. Setting up results

1. General settings

You can set the following settings: "Printing results when printing sample information, calibration and QC results" and "Dilution elements, expansion calculation results go beyond the linear range label".

2. X8 Settings

You can set the following settings: "Do not display concentration for test results with sampler errors" and "Display a label for results that goes beyond the measurement range".

In General settings or X8 settings, select the desired option and click <Save> to change the information.

5.7.1.6 Purpose of the analysis

Click <Purpose of the analysis> to go to the [Purpose of the analysis] screen, as shown in the figure below. Using the function, you can assign analyses to the appropriate tools.


Fig. 5.7-13. Purpose of the analysis

1. Purpose of the analysis

Select the Analysis module to which you want to assign the analysis, and select the analysis from the list of unassigned analyses. Click **<Add>** to assign the analysis to the tool. Click **<Add All>** if you want to assign all unassigned analyses to one tool at once.

2. Deleting the analysis

Select the Analysis module from which you want to delete the analysis. Select an analysis from the list of unassigned analyses and click **<Delete>**. Click **<Delete All>** if you want to delete all assigned analyses at once.

Note.

If no analyses are assigned to the tool, it cannot be used to perform sample testing and QC verification, but calibration can be performed.

5.7.1.7 Preset reagent volume

Click **<Preset reagent volume.>** to go to the screen **[Preset reagent volume]**, as shown in the figure below. The function is designed to adjust the approximate amount of reagent that is used per day.


Fig. 5.7-14. Preset reagents volume

1. Preset reagent volume

Select "Module" and "Analysis" in the [Reagent Usage Information per day] area, enter a value in the "Number of tests per day" field, select "Required" and click **<Add>** to set the amount of reagents. For analyses marked as mandatory, select the "Preliminary confirmation" option on the main page of the software, and the software will automatically determine whether the amount of reagent is sufficient.

2. Changing the preset reagent volume

In the list of preset reagent volumes, select an entry and change the information in the [Reagent Usage Information for the day] area on the right. Click **<Save>**.

3. Removal of the preset reagent volume

In the list of preset reagent volumes, select an entry and click **<Delete>** to delete it.

5.7.1.8 Warning for the reagent

Click **<Warning for the reagent>** to go to the screen [Warning for the reagent], as shown in the figure below. Using the warning settings for the reagent, you can specify whether the user will be shown a warning about the remaining tests for the reagent in the current analyzer. If the number of tests for the reagent is less than or equal to the specified amount, then in the area [Warning for the reagent] a message will appear on the main page stating that the reagent volume is insufficient.


Fig. 5.7-15. Warning for the reagent

Enter the number of remaining tests and click **<Save>** to specify the number of remaining tests at which the system displays a warning for the reagent.

5.7.2 System Settings

The second area on the [Settings] screen is the [System Settings] area, which contains the following settings: "Control. user", "Language", "Network connection", "Speakers", "Alarm limitation", "Module settings", "Version information" and "Backup copy/View". The following provides basic information about the various system settings.

5.7.2.1 Version Information

Click the **<Version Information>** button to go to the [Version Information] screen shown in the figure below. Information about the above-installed computer and the PLC version will be shown.


Fig. 5.7-16. Version Information

5.7.2.2 User Management

Click <User manag.> to go to the [User manag.] screen, shown in the figure below. The software includes user management functions that are used by authorized system administrators to assign rights to other users. The function of adding new users, changing the name, password and user rights is supported.


Fig. 5.7-17. User Management

1. Adding a user

In the [User Information] area, enter information about the new user and click the <Add> button.

2. Changing user data

In the [User List] area, select the user whose information you want to change, edit the information in the [User Information] area, and click the <Save> button to complete the data change.

3. Deleting users

In the [User List] area, select the user whose record you want to delete and click the <Delete> button to complete the deletion.

5.7.2.3 Language

Click the <Language> button to go to the [Language] screen shown in the picture below. After changing the language, the software will display the text in the appropriate language.


Fig. 5.7-18. Language

5.7.2.4 Network connection

To go to the screen [Network connection] (see the picture below) click the <Network connection> in the [System Settings] area.


Fig. 5.7-19. Network connection

The software supports two-way data exchange with the LIS server of the medical institution. Analysis information can be retrieved from the health care facility's LIS server by sample ID or barcode; analysis results are sent to the health care facility's LIS server.

1. Network Connection Settings

- Select the "No connection" option to work without using the functions available with a network connection.
- Select the COM port connection and the number of the corresponding COM port, as well as the transfer rate in the [COM Port Setup] section to complete the procedure.
- Select a TCP/IP connection and specify the appropriate IP address and port in the [Network Settings] section.

2. Other network connection settings

If necessary, the user can configure host verification, automatic operations, automatic uploading of results, and waiting time/number of resend attempts. When all settings are completed, click the **<Save>** button to complete the network connection settings procedure. After successful connection, the LIS status indicator located on the action bar at the bottom of the screen will turn white.

5.7.2.5 Speakers

Click the **<Speakers>** button to go to the **[Speakers]** screen shown in the figure below. This screen is used to enable/disable the sound signals through the speakers of the computer or analyzer in case of a warning or error, as well as to select the type of sound signal.


Fig. 5.7-20. Speakers

To go to the **[Signal Selection]** screen, click the **<Signal Selection>** button in the **[Speaker Setup]** area. Then select the audio file from the list of signals on the **[Signal Selection]** screen. To listen to an audio file, click the **<Play>** button, and to stop playback, click the **<Stop>** button. Click the **<OK>** button to complete the selection of the audio signal type.

Click the **<Save>** button on the **[Speakers]** screen to save the settings.

5.7.2.6 Limitation of the alarm signal

Click the **<Alarm Limit>** button to go to the **[Alarm Limit]** screen shown in the figure below. This screen is used to set up alarms that do not need to be set during the analysis.


Fig. 5.7-21. Limitation of the alarm signal

5.7.2.7 Module Settings

Click the **<Module Settings>** button to go to the **[Module Settings]** screen shown in the figure below. On this screen, you can set the settings for various modules.

When the SPS03 system is selected, you can specify the following: whether the module needs to perform sample analysis, scan the barcode of the sample, save the manually specified number, detect the tube, warn about the absence of a barcode on the tube, automatically download the analysis data for samples whose details are specified manually, and show the LC verification number in the sample registration area.

When the X8 system is selected, you can set the following in the general settings: showing information about how well the reagent is mixed; cooling in sleep mode; pump speed of the flushing unit.


Fig. 5.7-22. Module Settings

5.7.2.8 Backup copy/View

Click the **<Backup copy/View>** to go to the **[Backup copy/View] screen** shown in the figure below. The user can create a backup copy of the data or restore it.


Fig. 5.7-23. Backup copy/View

1. Exporting data

Select "End Date" and "Save Path" in the **[Create a backup copy and delete data]** area. Click **<Backup copy>** to export the database data and create a backup copy of it. Note. You can only export data for the last three days. You can also check the box "After creating a backup, all data with a date earlier than the end date will be deleted from the current database", and all exported data will be deleted from the current database.

2. Importing data

Select the "File Path" in the **Importing backup copy data]** area and click **<Import>**. The backup data will be imported into the database. After importing a backup copy of the data, you can select an archived database to view it on the **[Result]** screen, **[Res. of the calibr.]** or **[QC result]**.

5.7.3 System maintenance

The third area on the **[Settings]** screen is the **[System Maintenance]** area. It can be used to perform initialization and maintenance, as well as to troubleshoot components. The following provides basic information about the functions of the system maintenance area.

5.7.3.1 Initialization

Click the **<Initialization>** button to go to the **[Module Initialization]** screen. This module is used to initialize the analyzer, as shown in the figure below.


Fig. 5.7-24. Initialization

1. The device is fully equipped

Select [**Equipment.**] in the list on the left and click **<Initialization>** to initialize the device in full configuration.

2. Characteristics

Select the name of the device module from the list on the left and click the button on the right side to perform the appropriate actions, including those listed below.

Initialization of the tool.

Resetting the tool settings. When performing this operation, the cuvettes are not removed from the analyzer.

Sending a lower-level command. Enter the desired command in the [Low Level] field, then click the **<Send>** button.

Washing. Enter the number of flushes for the sampler, reagent dispenser holder, flushing unit and chamber. Click **<Flush>**. The tool will flush the relevant components.

Restarting the shaker. Allows you to restart the shaker of the tool after stopping.

5.7.3.2 Service

Press **<Maintenance>** to go to the [Maintenance of the device], where you can select a module and perform the appropriate maintenance action.

Select the chemiluminescence immunoassay module and click on the **<Daily>** tab to go to the daily maintenance tab. Follow the maintenance steps according to the suggested procedure. If it is necessary to perform a one-button cleaning, press **<Start>** for the appropriate procedure, after which the corresponding function will be performed.

The weekly and monthly maintenance screens have similar functions. Their detailed description is not provided in this document. For detailed information about daily, weekly and monthly maintenance, see Chapter 6, "System Maintenance".


Fig. 5.7-25. Daily maintenance of the chemiluminescence immunoassay module

5.7.3.3 Troubleshooting of components

Click <Troubleshooting of components> to go to the [Troubleshooting] screen. This screen also appears when the analyzer crashes or its components fail, as shown in the figure below.


Fig. 5.7-26. Troubleshooting of the chemiluminescence immunoassay module

1. Troubleshooting of the chemiluminescence immunoassay module

Click on the name of the tool on the left side of the dialog box to go to the troubleshooting screen of this tool. In the component initialization area above the dialog box, select one or more components to be separately initialized, and then click <Send> to start initializing the selected components.

If you want to pause the analysis, click the <Suspend cycle> button, and the analyzer will be interrupted.

Click the **<Stop>** button to stop all tests. Click the **<Continue>** button to continue performing analyses that can be performed after component initialization is complete.

Note.

Do not leave the system in a cycle suspended state for an extended period of time, as this may affect the results of some tests.

2. Elimination of failures of the sample processing system

Click on the name of the sample processing system on the left side of the dialog box to go to the Troubleshooting system screen and start the transportation device.

When will the sample processing system be in a testing state:

- Click **<Suspend Transport>**. The **<Suspend Transport>** button will be replaced by the **<Resume transport>**. The system will interrupt the output of the sample tripods from the sample area to the tool.
- After you click **<Resume transport>**, this button will change to the **<Suspend transport>**, and the system will continue to bring tripods with samples from the sample area to the tool.
- After you click **<Tripod Return>**, the tool will stop loading the samples and return the tripods to the sample area.

Note.

Here, the SPS03 (Sample Transfer System) includes a sample area and a sample conveyor for each module.


Fig. 5.7-27. Troubleshooting of components of the sample processing system

5.8 Description of [Home screen] screen functions

Click the <Home Screen> button on the main menu to go to the [Home Screen] menu and access its functions as shown in the figure below. The screen is divided into two sections: manual of operations and module monitoring. The following provides basic information about the components of the home screen.


Fig. 5.8-1. Home Screen

5.8.1 Manual of operations

The manual on the home screen shows the set reagent volume, reagent warning, utilities, calibration, maintenance status and settings for automatic sleep mode.

5.8.1.1 Preset reagents volume

As described in the chapter 5.7.1.7 "Preset reagent volume", you can check the "Preliminary confirmation" box after you have set the approximate volume of the reagent used and before the analyzer starts performing the analysis. The software will calculate the total volume of the reagent loaded into the analyzer and compare it with the preset volume.

Table. 5.8-1. Button colors <Preset reagent volume>

Colour	Description
Green	All preset reagent volumes are normal.
Yellow	The volume is preset for the reagent, and the current volume is less than this value.
Red	There is no reagent for which the volume is preset.

Press <Preset reagent volume> to go to the [Reminders of preset reagent volume] screen, where you can view information about a missing reagent or a reagent with a missing volume.


Fig. 5.8-2. Reminders of preset reagent volumes

5.8.1.2 Warning for the reagent

Reagent warnings are displayed for various reagent statuses, including when the reagent volume is insufficient to perform the remaining tests or the reagent's expiration date has expired.

Table. 5.8-2. Button Colors <Warning for the reagent>

Colour	Description
Green	Normal condition, all reagents are available.
Yellow	- The remaining reagent tests are in the range between >0 and ≤ 10. - The remaining reagent tests are more than zero and less than the amount specified for displaying the warning (or equal to it). For more information, see in the 5.7.1.8 "Warning for the reagent" section. - The shelf life of the reagent has expired.
Red	The number of remaining reagent tests is zero.

Click <Warning for the reagent> to go to the [Reagent: Warnings and reminders] screen, where you can view information about a missing reagent or a reagent with a missing volume.


Fig. 5.8-3. Warnings and reminders about reagents

5.8.1.3 Utility

The screen is used to display information about the system check of the chemiluminescence immunoassay module.

Table. 5.8-3. Button Colors <Utility>

Colour	Description
Green	The results of the system check are normal for the current day.
Yellow	The results of the system check for the current day are out of the norm.
Red	The system checks have not been completed on the current day.

Click <Utility> to go to the [Utility Reminder] screen, where you can view information about a specific system check alert.

Tap <Registration> on the [Utility Reminder] screen to go to the [Registration of data for system verification] screen. For more information, see in the section 5.6.1.1 "Registration of data for system verification").


Fig. 5.8-4 Utility Reminders

5.8.1.4 Calibration

The calibration status of all reagents loaded into the analyzer is displayed here.

Table. 5.8-4. Button Colors <Calibration>

Colour	Description
Green	All reagents have a valid calibration record or are in the process of calibration
Yellow	The recording of one of the reagents is new and unconfirmed.
Red	One of the reagents has not been calibrated, or its calibration has expired.

Tap <Calibration> to go to the [Immunoassay Calibration Reminder] screen, where you can view

information about a specific warning. On the **[Immunoassay Calibration Reminder]** screen, select reagent analysis from the list and click **<One-click Registration>** to en masse calibrations for this analysis.


Fig. 5.8-5. Calibration Reminders

5.8.1.5 Service

Allows you to view the maintenance status of all modules.

Table. 5.8-5. Button Colors <Maintenance>

Colour	Description
Green	The daily, weekly and monthly maintenance records for all modules have not expired, and the system does not remind you of the need for maintenance.
Yellow	The daily, weekly and monthly maintenance records for all modules have not expired, and the system reminds you about the maintenance.
Red	There are no records of daily, weekly, or monthly module maintenance, or a certain module has expired.

Click **<Maintenance>** to open the **[Maintenance Information]** dialog box to view detailed information.


Fig. 5.8-6. Maintenance Information

1. Setting the time

On the [Maintenance Information] screen, press <Set time.> to open the [Set Up Maintenance Reminders by Time] dialog box. For the "Daily Maintenance" module, select "Do not remind", "Remind when loaded" or "Remind at specified time" and click <Save>. Perform the same action for weekly and monthly maintenance.

2. Maintenance calendar

In the [Maintenance Information] dialog box, click <Calendar> to go to the [Maintenance Calendar] dialog box, where you can view information about daily maintenance.

3. Maintenance log

In the [Maintenance Information] dialog box, click <Log> to go to the [Maintenance Records] dialog box, where you can view the completed and ongoing maintenance work each day.

5.8.1.6 Automatic transition to sleep mode

The automatic sleep mode function allows you to set the settings for sleep mode and startup of the analyzer. Click <Automatic transition to sleep mode > to open the [Instrument Hibernation Settings] dialog box, where you can specify which components to hibernate, the time to hibernate, and the time to exit hibernation. Press the <Sleep mode> button to put the analyzer into sleep mode. Click <Start> and the analyzer will return to standby mode after initialization.

5.8.2 Module monitoring

The module monitoring function allows you to view the current status of each module. If the module enters a warning state, the background color of the module name in the monitoring area and the main menu changes to yellow. If the module enters a state where an error warning is displayed, the background color of the module name in the monitoring area and the main menu changes to red. If the module is in a normal state, the name of the module in the monitoring area and the main menu is shown on a green background. During the assessment of the tool status, a deviation from the normal state in the processing of the reagent and transaction will not change the background color of the module name in the monitoring area and the main menu. To view information about a specific status, click on the corresponding status button.

1. Status of the chemiluminescence immunoassay module

Click on the icon of the chemiluminescence immunoassay module to go to the main page of the module, where you can view the current status of the tool (see the section 5.9.1 "Status of the chemiluminescence immunoassay module"). In the status section of the chemiluminescence immunoassay module, you can find out the status of the following aspects: consumables, reagents, waste, environment, transaction processing and tools. The information is given in descending order of priority, including error information, warnings, and normal status messages. The three states with the highest priority are shown in the monitoring section of the module.

2. Status of the sample processing system

Click on the sample processing system icon to go to the home screen and view the current status of the system (see the section 5.9.2 "Status of the sample processing system"). The status of the sample processing system is shown in the monitoring section of the sample processing system as the status of the tool.

5.9 Description of the functions of the [Status] screen

Click on the icon of the corresponding tool on the main page (or select this tool from the main menu to go to the main page of the tool), as shown in the figure below. This screen displays the status of the analyzer components.

5.9.1 Status of the chemiluminescence immunoassay module


Fig. 5.9-1. Status of the chemiluminescence immunoassay module

The main page of the chemiluminescence immunoassay module shows the status and information about the processing of transactions for the environment, waste, consumables and parts of the tool components.

- If the component's health indicators are within the specified range or the component's condition is normal, the icon is displayed in green.
- If the component status indicators have reached a warning value, the icon is displayed in yellow.
- If the component's health indicators are outside the specified range or the component's condition is unsatisfactory, the icon is displayed in red.

To view specific values, click on the corresponding icon to go to the detailed information screen. This section discusses the functions associated with the states of the chemiluminescence immunoassay module.

5.9.1.1 Environmental parameters

In the environmental information area on the main tool screen, you can view information about temperature, humidity and voltage for different components of the tool. The following provides basic information about environmental parameters.


Fig. 5.9-2. Environment

Press the **<Temperature>** button on the home screen of the tool to go to the **[Temperature]** screen shown in the figure below. This screen is used to view the temperature or ambient temperature of each component of the analyzer in real time; temperature units can be adjusted. The value displayed on the icon corresponds to the actual value recorded in real time. If the indicator is within the normal range, the icon is displayed in green; if the indicator is outside the normal range, the icon is displayed in red. Screens with voltage and humidity indicators are arranged in the same way as the screen with temperature indicators.


Fig. 5.9-3. Temperature

5.9.1.2 Wastes

In the waste status information area of the tool's home screen, you can view the status of spent cuvettes, disposable tips, and fluid, as shown below. The following provides basic information about the waste conditions.


Fig. 5.9-4. Wastes

1. Used cuvettes and disposable tips

In the Waste information area, click **<Used Cuvettes>** or **<Disposed Tips>** to go to the **[Used Cuvettes and Tips]** screen, as shown in the figure below. You can view in real time the number and condition of used cuvettes and disposable tips, the total number of used disposable tips and empty tip containers, as well as information about the condition of empty cuvettes.


Fig. 5.9-5. Used cuvettes and tips

- If the number of used disposable tips is less than 250, the icon is displayed in green; if the number of used disposable tips is between 251 and 350, the icon is displayed in yellow; if the number of used disposable tips is more than 350, the icon is displayed in red. The user can set a value range at which a warning is displayed about the number of disposable tips used. For more information, contact the maintenance specialists.
- If the number of spent cuvettes is less than 2400, the icon is displayed in green; if the number of spent cuvettes is from 2401 to 3000, the icon is displayed in yellow; if the number of spent cuvettes is more than 3000, the icon is displayed in red.

Click on the icon  to change the number of relevant waste to zero.

Warning

1. After emptying the waste container, discard the amount of relevant used materials, otherwise the quantity information will not change.
2. Since the substances in the waste can lead to infection, they must be placed in a waste bag. Waste should be disposed of in accordance with local legal and regulatory requirements.

2. Liquid waste

In the Waste information area, click **<Liquid Waste>** to go to the **[Liquid Waste]** screen, as shown in the figure below. The function allows you to show the level status of liquids for high and low concentrations of liquid waste in real time.


Fig. 5.9-6. Liquid waste

If the liquid level is normal, the icon is displayed in green; if the liquid level meets the limit value, the icon is displayed in red.

5.9.1.3 Consumables

In the consumables status information area of the main tool screen, you can view information about the status of disposable tips, cuvettes, flushing liquid, starting reagents, cuvette return area and cooling water tank, as shown in the figure below. The following provides basic information about the states of consumables.

Fig. 5.9-7. Consumables

1. Disposable tip

In the Consumables status information area of the main tool screen, click **<Disposable Tip>** to go to the **[Disposable Tip]** screen, which is shown in the figure below. The screen that opens is used to display the current consumption status of disposable tips and tip blocks in real time, as well as to indicate the number of disposable tips for the current tip block.


Fig. 5.9-8. Disposable tips

- 1) The following is a description of the states of disposable tips shown on the home screen of the chemiluminescence immunoassay module.

If the number of disposable tips is >96, the icon is displayed in green; if the number of disposable tips is between 1 and 96, the icon is displayed in yellow; if the number of disposable tips is 0, the icon is displayed in red.

- 2) The following states of the disposable tip block are possible.

Status	Description
	Block of disposable tips is available
	Block of disposable tips is not available

- 3) The following conditions of disposable tips are possible.

Status	Description
	Disposable tips are present
	There are no disposable tips

- 4) Reset the number of disposable tips.

In the usage area of the current tip block located on the right side of the **[Disposable Tip]** screen, select the starting position for loading disposable tips and click . The current total number of disposable tips will change accordingly, indicating that the procedure for resetting the number of disposable tips has been completed.

Click on the icon  to empty all tool containers with tips.

Note.

1. The number of disposable tips after the reset must match the actual number of tips in the current disposable tip unit.
2. This function is only used if disposable tips are added manually or if there is no single disposable tip. If the initial loading position of disposable tips is not selected correctly, an error may occur when loading disposable tips.

2. Cuvettes

In the consumables status area of the main tool screen, click **<Cuvette>** or **<Cuvette Return>** to go to the **[Cuvette]** screen, which is shown in the figure below. This screen is used to display the number of cuvettes in the current cuvette block in real time, as well as the consumption of cuvette blocks in the storage area (in the sliding compartment for storing cuvettes and in the compartment for returning blocks). This screen is also used to indicate the number of cuvettes in the current cuvette block.


Fig. 5.9-9. Cuvette

- 1) Below is a description of the cuvette conditions, which are displayed on the main page of the chemiluminescence immunoassay module.

When the number of cuvettes is >364, the icon is displayed in green; when the number of cuvettes is between 1 and 364, the icon is displayed in yellow; when the number of cuvettes is 0, the icon is displayed in red.

- 2) The following storage area states are possible:

Status	Description
1	Cuvette blocks are available
2	There are no cuvette blocks

- 3) The following conditions of the compartment for used blocks are possible:

Status	Description
1	Cuvette blocks are available
2	There are no cuvette blocks

4) The following conditions of the cuvette are possible:

Status	Description
	The cuvettes are available
	There are no cuvettes

5) Reset the number of cuvettes.

In the information area about using the current cuvette, located on the right side of the [Cuvette] screen, select the initial position of the cuvette capture and press the button . The current total number of cuvettes will change accordingly, and the number of cuvettes will reset.

Click on the icon  to empty all the tool containers with cuvettes.

Note.

1. The number of cuvettes after the reset must correspond to the actual number of cuvettes in the current block.
2. This function is used only if the cuvette is added manually or if there is no single cuvette. If the initial position of the cuvette capture is selected incorrectly, an error may occur when capturing the cuvette.

3. Washing liquid and starting reagent

In the consumables status information area of the main tool screen, click <Washing liquid> to go to the [Washing liquid] screen, which is shown in the figure below. The remaining volume of the washing liquid is displayed on this screen in real time.


Fig. 5.9-10. Washing liquid

When setting the automatic mode of the tank for the washing liquid:

- When the washing liquid/clean water tank is $\geq 10\%$ of the full load, the icon will be displayed in green;

- When the washing liquid/clean water tank is less than 10% of the full load, the icon is displayed in yellow;
- When the fixed level in the tank is low, the icon is displayed in red.

When setting the manual mode of the washing liquid container:

- If the total volume is $\geq 20\%$ of the possible download volume, the icon is displayed in green.
- If the total volume is $< 20\%$ of the possible loading volume, and when the liquid level in the fixed tank is above the lower level, the icon is displayed in yellow.
- When the level of the flushing liquid in the stationary tank is low, the icon is displayed in red.

In the Consumables status information area of the main tool screen, click **<Starting reagent>** to go to the **[Starting reagent]** screen, which is shown in the figure below. This screen is used to view the total volume of starting reagents 1 and 2 in real time.


Fig. 5.9-11. Starting reagent

- If the remaining volume of the starting reagent 1 or 2 is $\geq 20\%$ of the possible loading volume, the icon of the corresponding reagent is displayed in green.
- If the remaining volume of the starting reagent 1 or 2 is $< 20\%$ and $\geq 4\%$ of the possible loading volume, the icon of the corresponding reagent is displayed in yellow.
- If the remaining volume of the starting reagent 1 or 2 is zero, the icon of the corresponding reagent is displayed in red.

5.9.1.4 Device status

You can view the status of the shaker and the disc of the tool in the "Status" area on the home screen of the tool, as shown in the figure below. The following provides basic information about the states of the tool.


Fig. 5.9-12. Tool States

1. Disk status

In the Consumables status information area of the main tool screen, click **<Disk Status>** to go to the **[Disk Status]** screen, which is shown in the figure below. With this function, you can monitor the cuvette positions and statuses on each panel component in various tools in real time.

In the central part of the **[Disk Status]** screen, the following components are displayed (in order from left to right and from top to bottom): the flushing unit, the reaction disk, the chamber and the sample disk. The current states of the cuvettes on the disk are indicated by different colors. There is a list of symbols on the left side of the screen. The list of messages for the current analysis performed for the cuvette at the selected position is displayed on the right side of the screen.


Fig. 5.9-13. Disk status

2. Shaker status

In the tool status information area on the main tool screen, you can monitor the status of the left and right shakers of each tool. If the speed of the shaker is within the normal range, the icon is displayed in green. If the speed of the shaker is outside the acceptable range, the icon is displayed in red.

5.9.1.5 Processing operations

In the Transaction status information area on the main tool screen, you can view information about the upcoming calibration, the last BGW check, and the last maintenance of each tool, as shown in the figure below. The following is a summary of the basic information about the transaction processing states.


Fig. 5.9-14. Processing operations

1. Awaiting calibration

In the Transaction status information area on the main tool screen, click **<Pending Calibration>** to go to the **[Reagent]** screen of the corresponding tool. The function is designed to calibrate the reagent and perform other actions.

2. Last check of BGW

In the transaction status information area on the home screen of the tool, click **<Last check of BGW>** to go to the **[System checks]** screen. The function is designed to perform data logging for system verification and other actions with the tool.

3. The time of the last service

In the consumables status information area of the main tool screen, click **<Last service>** to go to the **[Maintenance Records]** screen, which is shown in the figure below. The function allows you to view, delete, export and print tool maintenance records.


Fig. 5.9-15. Maintenance records

5.9.2 Status of the sample processing system

Click on the sample processing system icon on the main page (or select this icon in the main menu to go to the main page) to go to the main page of this system. Press **<SPS Status>** to go to the **[Sample Processing System Status]** screen. You can view the status of the sample stand in the sample area and on the tripod.

If the tool fails to add samples during analysis, and the tripod with samples remains in the same

loading position of the tool, then you can double-click on the sample tripod icon and press <Unlock> to unlock the tripod and return it to the sample area.

The color of the left sample area indicator in the software corresponds to the sample area indicator of the tool. The values of the indicator colors are described in the section 4.2.7.


Fig. 5.9-16. Status of the sample processing system

6.1.2 Materials to be prepared

Clean cotton cloth: to clean each component.

75 percent alcohol: for cleaning various parts.

Cotton swab: used to clean the reagent dispenser.

Pliers (or scissors): these can be used to cut the clamps on the tubes.

Clamps: for tightening the tubes.

6.1.3 Cuvettes

Cuvettes are used to transfer samples and reagents needed for the reaction. All cuvettes should be purchased only from Shenzhen New Industries Biomedical Engineering Co., Ltd. The use of cuvettes from other manufacturers is prohibited. Our company is not responsible for inaccurate results due to the use of cuvettes manufactured by other companies.

6.1.4 Disposable tip

Disposable tips are used to transfer or temporarily store samples during analyzer assays. Only disposable tips purchased from Shenzhen New Industry Biomedical Engineering Co., Ltd. may be used. If other disposable tips are used, the accuracy and reproducibility of the assay results may be insufficient. SNIBE accepts no responsibility for inaccurate results when using disposable tips not purchased from SNIBE.

6.1.5 Washing liquid

Washing liquid is used to clean all components of the analyzer. Only washing liquids purchased from Shenzhen New Industry Biomedical Engineering Co., Ltd. may be used. If other flushing fluids are used, the cuvettes, reagent pipettes, and tubing may not be cleaned sufficiently, and the assay mixture in the cuvette may not be cleaned sufficiently, adversely affecting the accuracy and reproducibility of the assay results. SNIBE assumes no responsibility for inaccurate results when cleaning cuvettes using washing liquid not purchased from SNIBE.

6.1.6 Starting reagent

The starting reagent can react with the luminescence tag in the supported chemiluminescence reagent. The reaction is accompanied by a light signal of a certain intensity required for analysis by the analyzer. Only starter reagents purchased from Shenzhen New Industries Biomedical Engineering Co., Ltd. may be used. If other starting reagents are used, the accuracy and reproducibility of the assay results may be insufficient. SNIBE accepts no responsibility for inaccurate results when using starting reagents not purchased from SNIBE.

6.1.7 Solution for cleaning system tubes

The basic composition of the system tube cleaning solution can dissolve protein, which is required to maintain the analyzer in proper condition. Only the system tube cleaning solution purchased from Shenzhen New Industries Biomedical Engineering Co., Ltd. may be used. If another liquid is used to flush the tubing, the effectiveness of the tubing flushing may be reduced, which may result in poor accuracy and reproducibility of the assay results. SNIBE assumes no responsibility for inaccurate results when using system tube cleaning solution not purchased from SNIBE.

6.1.8 Purified water

For daily operation and maintenance of the analyzer, purified water with a conductivity of 1 $\mu\text{S}/\text{cm}$ or less should be used. When using the purified water unit, regular maintenance and inspection should also be ensured, please refer to the purified water instruction manual for details or contact the manufacturer or retailer of the purified water unit.

6.2 Daily maintenance

Materials needed:

Clean waste bag, clean cotton cloth

Maintenance Procedures:

1. Use a clean cotton cloth to clean the surface of the analyzer.
2. Empty the two waste bags of used cuvettes and disposable tips and replace the bags with new ones. Empty the liquid waste tank.
3. Ensure that the volume of starting reagent is sufficient for the next day's assays, add or replace the starting reagent as necessary.
4. Ensure that the volume of washing liquid is sufficient for the next day's assays, when manual preparation mode is selected, then add washing liquid if necessary.
5. Perform one-button washing once.

6.3 Weekly maintenance

Materials needed:

Cotton swabs, alcohol.

Maintenance Procedures:

1. Using an alcohol soaked cotton swab, clean the tip of the sampler holder needle in the fluid handling module, the outside of the reagent probe, the outside of the rinse needle, and the liquid waste needle in the rinse assembly and chamber.
2. Clean the surface of the cuvette and disposable tip storage compartments with a clean cotton cloth.
3. Perform one-button washing twice.

6.4 Monthly maintenance

Materials needed:

Sodium hypochlorite solution with effective chlorine content of 0.05%~0.1% (i.e. 500mg/L~1000mg/L), clean water, system tube cleaning solution, light signal test solution.

Maintenance Procedures:

1. Replace the starting reagent with 40□ clean water and rinse the chamber 30 times, then replace the clean water with the starting reagent and rinse the chamber another 30 times.
2. Use pipe flushing fluid and perform a one-button washing operation once.
3. Empty the washing liquid tank and the liquid waste tank.
 - 1) Empty the washing liquid tank and the liquid waste tank.
 - 2) Add 2 liters of sodium hypochlorite solution with an effective chlorine content of 0.05%~0.1% (i.e., 500 mg/L~1000 mg/L) to the barrel, cover with a lid, shake back and forth, and leave for 30 min;
 - 3) Pour the cleaned sodium hypochlorite solution and add clean water (2 liters each time) to the cleaning bucket at least 3 times.
4. Fill the bucket with cleaning solution to spray the pipette dispenser and flushing assembly, 30 times each;
5. Place the container of light test liquid in the upper right corner of the sample compartment. Perform the BGW test and the light signal test once each. Review the results and record them.
6. System verification requirements and results:
 - 1) BGW 100–1200 coefficient of variation ≤ 10 %
 - 2) LC450 000–650 000 coefficient of variation ≤ 3 %
7. For specific expected values, see the instructions for use of the LC test reagent.

6.5 Prolonged non-use

If the analyzer will not be used for more than 5 days, perform the following steps:

Replace starting reagents 1 and 2 with clean water.

Replace the washing liquid solution in the washing liquid reservoir with clean water.

Select the appropriate chemiluminescence immunoassay module. from the [Module Initialization] screen, press the <Rinse> button, and flush the sampler, reagent dispenser holder, wash assembly, and chamber 30 times, filling the tubes with clean water.

When you need to use the analyzer again, follow these steps:

1. First remove the bottle of starting reagent and the container of the washing liquid. Ensure that the tube openings are not contaminated. Click the [Module Initialization] button to select the appropriate chemiluminescence immunoassay module and click the [Washing] button on the screen and flush the sampler, reagent dispenser holder, flushing assembly, and chamber 30 times, then empty the tubes.

2. Add starting reagents 1 and 2.
3. Add washing liquid.
4. Press the [Module Initialization] button to select the appropriate chemiluminescence immunoassay module and click the <Start> button on the [Module Initialization] screen and wash the sampler, reagent dispenser holder, wash assembly, and chamber 30 times each. Before starting the analysis, make sure that the results of the system check meet the requirements.

Chapter 7 The main malfunctions and ways to eliminate them

7.1 Classification of faults

The faults can be divided into two main categories: alarmed and non-alarmed. Depending on the severity of the malfunction, a warning alarm, a sample stop loading alarm, and an analyzer suspend alarm are issued. The following table lists the actions performed by the analyzer in the event of an alarm.

Table. 0-1. Analyzer responses to an alarm

Alarm type	Analyzer action
Warning alarm	[Immunity Module] Sample loading and measurements will not be paused, the analyzer will continue to run the assays
	[Sample Area] The sample area continues to operate.
Alarm with sample loading stopped	The introduction of the new sample will be stopped, but samples already added to the cuvettes will continue to be analyzed
Alarm with shutdown	All analyzer operations will stop and the analysis results will be unreliable

7.2 Faults not accompanied by an alarm and how to rectify them

Some deviating events may not be accompanied by an alarm. Such events can be detected by checking control procedures, reviewing results, and checking analyzer status.

Table. 7.2-1. Data deviations not accompanied by an alarm and how to correct them

Signs of malfunction	Possible reason	What needs to be checked	Decision
Abnormal background check result: CV or RLU values outside the permissible range	Either the washing liquid has been stored too long or the problem with the water quality in the system, or a configuration error led to a malfunction related with washing liquid	Ensure that the washing liquid is not too viscous, the washing liquid container is free of yellow colored contaminants and flakes, the quick-switch port is clean, and the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$.	Clean the washing liquid container according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{cm}/\text{cm}$ and re-prepare the washing liquid according to instructions. Rinse the washing unit 30 times.
	The needle of the washing unit is dirty or there was a malfunction during the injection of the washing liquid	Make sure that the needle of the washing unit is clean, the four needles of the washing unit supply a constant amount of washing liquid, that the supply tube is not bent, clogged or leaking, and the locking screws of the solenoid valve and the piston pump are tightened.	Wipe the needle of the washing unit or contact technical support
	Malfunction of the peristaltic liquid waste pump of the washing unit, blockage or contamination of the liquid waste removal needle of the washing unit, the position of the washing liquid pump is incorrect, the needle tube is not elastic	Check the liquid flow rate in the liquid waste tube, the presence of residual liquid in the cuvettes after pumping, and contamination of the liquid waste disposal needle.	Wipe the needle of the washing unit or contact technical support.

	The flushing cuvette of the liquid waste removal pump needle is dirty, the tube is clogged, or the position of the washing angle is incorrect	Check for contamination of the washing cuvette of the liquid waste removal pump needle and the drain rate of the washing liquid. Make sure that the flushing position of the liquid waste removal needle in the washing cuvette is correct.	Contact technical support
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support
	An unacceptable amount of starter reagent was introduced by the starter reagent pump	Ensure that the starting reagent is sufficient, the pump does not produce unusual sounds and the starting reagent tube is not displaced or bent	Contact technical support
	Pollution of the cuvette	Check how long the cuvettes have not been used and ensure that they are free of dust and other contaminants. Check the environment of use or storage for contaminants.	Replace the cuvettes and clean the cuvette and tip storage compartments
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	The analyzer was used for a long time, and tube maintenance was not performed properly	Check the tubes for signs of wear and dirt	Run the tube flushing function once. When performing this procedure, it is necessary to use a flushing liquid for the tube system.
	The shelf life of the starting reagent has expired	Check the shelf life of the starting reagent and make sure that the storage conditions of the starter reagent meet the necessary requirements	Replace the starting reagent and flush the flushing assembly 30 times

Signs of malfunction	Possible reason	What needs to be checked	Decision
	Leakage in the tube of the reagent dispenser holder	Make sure that after cleaning the reagent dispenser, there are no air bubbles or water droplets at its tip, there are no air bubbles in the tube of the reagent dispenser holder, and also check the connection points of the tubes (solenoid valve, piston pump, etc.) for the presence of crystallization that could form due to leakage	Check the locking rings at the valve junction for leaks, make sure that there is enough flushing liquid and the reagent dispenser holder is in the correct position for flushing. Then select the reagent dispenser holder flushing option on the "Washing" screen to perform the flushing. If the leak continues, contact technical support
	Failure of the piston pump of the reagent dispenser holder, valve or flushing pump	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the flushing pump is not damaged and the liquid is not sprayed	Contact technical support
	Liquid level detection error, LC test solution not added	Check if there are LC test results close to 0	Contact technical support
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support
	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	Disposable tips or cuvettes are dirty	Check how long it has been since cuvettes or disposable tips have been used, and make sure there is no dust or other contaminants	Replace disposable tips or cuvettes and clean the cuvette and tip storage compartments

	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	Air leaks of the sampler holder pipeline	Check whether there are no air bubbles in the pipeline of the sampler holder, as well as whether there are leaks and crystallization of liquid in the pipe connection positions (solenoid valve, plunger pump, etc.).	Check if there is an air leak in the locking ring at the valve connection, check if the cleaning solution is not used up, check the correct position for cleaning the sampler holder, and then operate the interface for perfusion of the sampler holder. If the leak still exists, contact technical support.
	Malfunction of the plunger pumps of the sampler holder and the cleaning pumps or valves	Check the alarm recording to make sure that the plunger pump can be initialized, and check if the inner membrane of the cleaning pump is damaged and if liquid is sprayed.	Contact technical support.
	The analyzer was used for a long time, and tube maintenance was not performed properly	Check the tubes for signs of wear and dirt	Run the tube flushing function once. When performing this procedure, it is necessary to use a flushing liquid for the tube system.

Signs of malfunction	Possible reason	What needs to be checked	Decision
High coefficient of variation of the analysis results	Either the washing liquid has been stored too long or the problem with water quality in the system, either a configuration error resulted to a malfunction related with washing liquid	Ensure that the washing liquid is not too viscous, the washing liquid reservoir is free of yellow colored contaminants and flakes, the quick-switch port is clean, and the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$.	Clean the washing liquid container according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{cm}/\text{cm}$ and re-prepare the solution of the washing liquid in the working concentration according to instructions. Rinse the washing unit 30 times.
	Leakage in the tube of the reagent dispenser holder	Make sure that after cleaning the reagent dispenser, there are no air bubbles or water droplets at its tip, there are no air bubbles in the tube of the reagent dispenser holder, and also check the connection points of the tubes (solenoid valve, piston pump, etc.) for the presence of crystallization that could form due to leakage	Check the locking rings at the valve junction for leaks, make sure that there is enough flushing liquid and the reagent dispenser holder is in the correct position for flushing. Then select the reagent dispenser holder flushing option on the "Washing" screen to perform the flushing. If the leak continues, contact technical support
	Failure of the piston pump of the reagent dispenser holder, valve or flushing pump	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the flushing pump is not damaged and the liquid is not sprayed	Contact technical support

	Liquid level detection error, sample or reagent not added	Check your most recent test results (as many tests as possible) to see if your the RLU values are too high or too low	Contact technical support
	The needle of the washing unit is dirty or there was a malfunction during the injection of the washing liquid	Make sure that the needle of the washing unit is clean, the four needles of the washing unit supply a constant amount of washing liquid, that the supply tube is not bent, clogged or leaking, and the locking screws of the solenoid valve and the piston pump are tightened.	Wipe the needle of the washing unit or contact technical support
	Malfunction of the peristaltic liquid waste pump of the washing unit, blockage or contamination of the liquid waste removal needle of the washing unit, the position of the washing liquid pump is incorrect, the needle tube is not elastic	Check the liquid flow rate in the liquid waste tube, the presence of residual liquid in the cuvettes after pumping, and contamination of the liquid waste disposal needle.	Wipe the needle of the washing unit or contact technical support.
	The flushing cuvette of the liquid waste removal pump needle is dirty, the tube is clogged, or the position of the washing angle is incorrect	Check for contamination of the washing cuvette of the liquid waste removal pump needle and the drain rate of the washing liquid. Make sure that the flushing position of the liquid waste removal needle in the washing cuvette is correct.	Contact technical support
	The chamber or the light-permeable window is dirty	Check for liquid or crystals	Contact technical support

	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	Disposable tips or cuvettes are dirty	Check how long it has been since cuvettes or disposable tips have been used, and make sure there is no dust or other contaminants	Replace disposable tips or cuvettes and clean the cuvette and tip storage compartments
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
	Air leaks of the sampler holder pipeline	Check whether there are no air bubbles in the pipeline of the sampler holder, as well as whether there are leaks and crystallization of liquid in the pipe connection positions (solenoid valve, plunger pump, etc.).	Check if there is an air leak in the locking ring at the valve connection, check if the cleaning solution is not used up, check the correct position for cleaning the sampler holder, and then operate the interface for perfusion of the sampler holder. If the leak still exists, contact technical support.
	Malfunction of the plunger pumps of the sampler holder and the cleaning pumps or valves	Check the alarm recording to make sure that the plunger pump can be initialized, and check if the inner membrane of the cleaning pump is damaged and if liquid is sprayed.	Contact technical support.
	The analyzer was used for a long time, and tube maintenance was not performed properly	Check the tubes for signs of wear and dirt	Run the tube flushing function once. When performing this procedure, it is necessary to use a flushing liquid for the tube system.

	Air bubbles or foreign particles were detected in the sample or reagent	External inspection	Remove the air bubbles
	Insufficient centrifugation of serum	External inspection	Standardize the time and speed of serum centrifugation
	Test tube with the sample for analysis it is not installed in a tripod for samples or is not in an upright position	External inspection	Install the test tube with the sample for analysis in the sample stand properly

Signs of malfunction	Possible reason	What needs to be checked	Decision
Cross contamination	Either the washing liquid has been stored too long or the problem with water quality in the system, either a configuration error resulted to a malfunction related with washing liquid	Ensure that the washing liquid is not too viscous, the washing liquid reservoir is free of yellow colored contaminants and flakes, the quick-switch port is clean, and the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$.	Clean the washing liquid container according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{cm}/\text{cm}$ and re-prepare the washing liquid according to instructions.
	The needle of the washing unit is dirty	Check the needle of the washing unit for contamination.	Wipe the needle of the washing unit or contact technical support
	The flushing cuvette of the liquid waste removal pump needle is dirty, the tube is clogged, or the position of the washing angle is incorrect	Check for contamination of the washing cuvette of the liquid waste removal pump needle and the drain rate of the washing liquid. Make sure that the flushing position of the liquid waste removal needle in the washing cuvette is correct.	Contact technical support
	The reagent dispenser is dirty or damaged	Check the outer surface of the dispenser for crystals or serious damage to the coating	Wipe the reagent dispenser. If the dispenser is damaged, contact technical support for replacement
	Failure of the piston pump of the reagent dispenser holder, valve or flushing pump	Check if the piston pump initialization failure error is logged. Make sure that the inner membrane of the flushing pump is not damaged and the liquid is not sprayed	Contact technical support

Signs of malfunction	Possible reason	What needs to be checked	Decision
The RLU value too low according to the analysis results (or even equal to 0)	Liquid level detection error, sample or reagent not added	Check your most recent test results (as many tests as possible) to see if your the RLU values are too high or too low	Contact technical support
	Insufficient sample quantity	External inspection	Add the sample
	Atypical supply of the starting reagent by the pump	Make sure that the starting reagent is sufficient, the pump does not produce unusual sounds, the starting reagent supply tube is not displaced or bent, and there are no crystals on its outlet	Contact technical support
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support
The RLU value too high according to the analysis results	The concentration of the sample exceeds the measurement range	—	Perform the dilution
	Either the washing liquid has been stored too long or the problem with water quality in the system, either a configuration error resulted to a malfunction related with washing liquid	Ensure that the washing liquid is not too viscous, the washing liquid reservoir is free of yellow colored contaminants and flakes, the quick-switch port is clean, and the electrical conductivity of the clean water used to prepare the washing liquid is less than 1 $\mu\text{S}/\text{cm}$.	Clean the washing liquid container according to the monthly maintenance procedure, use clean water with an electrical conductivity of less than 1 $\mu\text{cm}/\text{cm}$ and re-prepare the washing liquid according to instructions.

	The needle of the washing unit is dirty or there was a malfunction during the injection of the washing liquid	Make sure that the needle of the washing unit is clean, the four needles of the washing unit supply a constant amount of washing liquid, that the supply tube is not bent, clogged or leaking, and the locking screws of the solenoid valve and the piston pump are tightened.	Wipe the needle of the washing unit or contact technical support
	Malfunction of the peristaltic liquid waste pump of the washing unit, blockage or contamination of the liquid waste removal needle of the washing unit, the position of the washing liquid pump is incorrect, the needle tube is not elastic	Check the liquid flow rate in the liquid waste tube, the presence of residual liquid in the cuvettes after pumping, and contamination of the liquid waste disposal needle.	Wipe the needle of the washing unit or contact technical support.
	The flushing cuvette of the liquid waste removal pump needle is dirty, the tube is clogged, or the position of the washing angle is incorrect	Check for contamination of the washing cuvette of the liquid waste removal pump needle and the drain rate of the washing liquid. Make sure that the flushing position of the liquid waste removal needle in the washing cuvette is correct.	Contact technical support
	Disposable tips or cuvettes are dirty	Check how long it has been since cuvettes or disposable tips have been used, and make sure there is no dust or other contaminants	Replace disposable tips or cuvettes and clean the cuvette and tip storage compartments
	Light scattering from the chamber	Check whether the RLU values are abnormally high within a single series of assay results	Contact technical support
	PMT malfunction	Ensure that the dark current noise level is within normal limits	Contact technical support

	Insufficient centrifugation of serum	External inspection	Standardize the time and speed of serum centrifugation
--	--------------------------------------	---------------------	--

7.3 Faults not accompanied by an alarm and how to rectify them

If an alarm occurs, check and acknowledge the alarm in the analyzer, then follow the instructions in the dialog boxes that appear. If the error cannot be cleared or if another problem occurs, contact a service technician or authorized distributor.

Chapter 8 Environmental requirements

The analyzer must be manufactured and operated in an environmentally safe manner and in compliance with the rules and requirements of the applicable legislation concerning environmental protection requirements. At the end of its useful life, the analyzer must be disposed of in accordance with the instructions governing the disposal of this type of equipment. the analyzer must be disposed of in accordance with local biosafety regulations. the analyzer cannot be disposed of with household waste. The analyzer must be delivered to a specialized, accredited (licensed) collection point for the recycling of electrical and electronic equipment.

8.1 Disposal and destruction

The service life of the analyzer is 10 years. At the end of its useful life, the analyzer must be disposed of in accordance with the instructions governing the disposal of this type of equipment. If you have any questions regarding disposal, please contact the manufacturer.

Dispose of all materials in a safe acceptable manner that complies with SanPiN 2.1.3684-21 and SP 1.3.2322-08. Equipment disposal is carried out by licensed organizations licensed for this type of activity.

Certain substances in the analyzer must be disposed of in accordance with environmental regulations. Dispose of in accordance with local regulatory requirements.

All departments using analyzers must dispose of liquid waste in accordance with toxicity standards and the manufacturer's or authorized dealer's recommendations.

Waste fluids from patients with infectious diseases should be disposed of using the devices designated for this purpose.

Chapter 9 List of international standards applicable to MD

Point	Standard reference	Title
-------	--------------------	-------

1	EN ISO 14971	Medical devices Application of risk management to medical devices.
2	EN 62366-1	Medical devices Part 1 Design of medical devices with regard to usability
3	IEC/TR 62366-2	Medical devices Part 2 Guidance on the application of usability design to medical devices
4	EN ISO 13485	Medical devices Quality management systems Requirements for regulatory purposes
5	EN ISO 18113-1	In-vitro diagnostic medical devices, information provided by the manufacturer (labeling), part 1 terms, definitions and general requirements
6	EN ISO 18113-3	In-vitro diagnostic medical devices, information provided by the manufacturer (labeling), part 3 in-vitro diagnostic instruments for professional use
7	EN 62304	Medical devices. Software Life-cycle processes
8	EN ISO 15223-1	Medical Devices Symbols used in the marking of medical devices, on labels and in accompanying documentation Part 1 basic requirements
9	EN 13612	Evaluating the performance of the medical devices for in vitro diagnostics
10	EN 61010-1	Safety requirements for electrical equipment for measurement, control and laboratory use. Part 1. General requirements
11	EN 61010-2-101	Safety requirements for electrical equipment for measurement, control and laboratory use Part 2-101 Particular requirements for in vitro diagnostic medical equipment
12	EN IEC 61010-2-081	Safety requirements for electrical equipment for measurement, control, and laboratory use-Part 2-081: Particular requirements for automatic and semi-automatic laboratory equipment for analysis and other purposes
13	EN 61326-1	Electrical equipment for measurement, control and laboratory use ems requirements part 1: general requirements
14	EN 61326-2-6	Electrical equipment for measurement, control and laboratory use. Electromagnetic compatibility requirements. Part 2-6. Specific requirements for in vitro diagnostic (ivd) medical devices.
15	EN IEC 63000	Technical documentation for the evaluation of electrical and electronic products with regard to the restriction of hazardous substances

Chapter 10 Storage, transportation, packaging

This medical device is transported by all modes of transport in covered vehicles in accordance with the requirements and rules of transportation of goods applicable to each type of transport.

Disconnect the power cord from the mains before storing. Protect the Analyzer from dust and other contaminants during storage. It is recommended to store the analyzer in its original packaging. The analyzer must be individually packed before transportation.

Storage/transportation requirements

- During transportation, the analyzer must be kept in a strictly upright position.
- During transportation, the analyzer must be protected from moisture, water, strong vibrations and compression. Use caution when loading and unloading the unit.
- The analyzer should be stored in a well-ventilated room protected from chemicals, corrosive gases, and strong sunlight.
- Storage/transportation temperature -20~55 °C
- Relative humidity during storage/transportation $\leq 93\%$ RH
- Atmospheric pressure during transportation and storage 500 hPa ~ 1060 hPa

Appendix A. Operational characteristics

1.1. Accuracy and repeatability of dosing

Check the sample pipetting volume of 10 μl and 280 μl and the reagent pipetting volume of 10 μl and 450 μl , respectively. The results should correspond to the values in Table 1.

Table 1. Requirements for accuracy and repeatability of dispensing

Nominal dispensing volume (volume)/ μl	Requirements	
	Offset	Coefficients of variation (CV) (%)
10	Within $\pm 1 \mu\text{l}$	≤ 5
280	Within $\pm 5\%$	≤ 2
450	Within $\pm 5\%$	≤ 2

1.2 Accuracy and deviations of the reaction disk temperature

The temperature error of the reaction disk shall be $\pm 0.3 \text{ }^{\circ}\text{C}$ from the setpoint; deviations shall not exceed $0.1 \text{ }^{\circ}\text{C}$.

1.3 LC check components

1.3.1. Analyzer noise level

After switching on the analyzer and its transition to the stable operation mode, the noise level should be in the range of 100-1200, the coefficient of variation should be $\leq 10 \%$.

1.3.2. Linearity of the RLU measurement

If the order of magnitude of the RLU using the luminescence method is at least 3, the linear correlation coefficient (r) must be greater than 0.99 (or must be equal to this value).

1.3.3. Repeatability of the RLU measurement results

When using the luminescent method, the coefficient of variation (CV, %) should not exceed 5 %.

1.3.4. Stability of the RLU measurement results

When using the luminescence method, the variation of the RLU values shall not exceed $\pm 10 \%$.

1.4 Transfer

The transfer index should be $\leq 10^{-7}$.

1.5 Reproducibility of clinical assays within a single batch

The coefficient of variation (CV, %) should not exceed 5%.

Appendix B. Network interfaces

B.1 Protocol

The software utilizes ASTM E1394-97 and HL7 data transfer protocols.

All data transferred between the software and the LIS system is stored in the following path:
online.log.

B.2 Protocol ASTM

B.2.1 Description of control codes

Table B.2-1. Characteristics of the control code

Code	Description	Corresponding ASCII code
<ENQ>	Query	0x05
<ACK>	Confirmation of response	0x06
<STX>	Start of the document	0x02
<ETX>	End of the document	0x03
<CR>	Start of the line	0x0D
<EOT>	Completing data transfer	0x04

B.2.2 Description of the standard data transfer format

<ENQ><STX>content<STX><EOT>

For example:

<ENQ><STX>TEST<STX><EOT>

Table B.2-2. Characteristics of the standard data transmission format

Code	Description	Corresponding ASCII code
<ENQ>	Query	0x05
<STX>	Start of the document	0x02
T	The letter T	0x54
E	The letter E	0x45
S	The letter S	0x53
T	The letter T	0x54
<ETX>	End of the document	0x03
<EOT>	Completing data transfer	0x04

B.2.3 Description of the separators

The ASTM E1394-97 protocol provides 5 transmission data separators. For more information, see section 6.4 of the E1394 protocol. The value of each separator is shown in the table below.

Table B.2-3. Description of the separators

Separators	Description	Corresponding ASCII code
<CR>	Record separator	0x0D
	Field separator	0x7C
\	Repeat separator	0x5C
^	Component separator	0x5E
&	Special separator	0x26

B.2.4 Description of the message types

Table B.2-4. Types of the messages

Message label	Description
H	Recording the message header
P	Recording patient information
O	Recording analysis information
R	Recording the result
Q	Recording the query for information
L	Finalizing the message

B.2.4.1 Recording the message header (H)

Note.

The information in this section is provided in accordance with Section 7 of ASTM Protocol E1394-97. The message header entries must be at the beginning of all transmission entries. They contain basic information about protocols and basic data transmission.

Example:

H|^&||PSWD|MAGLUMI X8|||||LIS||P|E1394-97|20180717<CR>

Table B.2-5. Recording the message header

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
-------	-------	---------------------------------------	---------	----------------	-------------

1	7.1.1	Message type	H	1	Yes
2	7.1.2	Separator definition	\ ^ &	4	Yes
4	7.1.4	Password	PSWD	20	No
5	7.1.5	Sender name	MAGLUMI X8	20	Yes
10	7.1.10	Receiver name	LIS	20	Yes
12	7.1.14	Processing method	P	1	Yes
13	7.1.13	Protocol version number	E1394-97	10	Yes
14	7.1.14	Date	YYYYMMDD	14	Yes

Note.

- 1) The total number of fields in the protocol is 14. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software. The "|" separator should also be used to separate different fields if some of these fields do not require values.
- 2) The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

B.2.4.2 Recording of the patient information (P)

Note.

The information in this section is provided in accordance with Section 8 of ASTM Protocol E1394-97. Information related to the patient should be presented using this type of record.

Example:

P|1|101|||ABC|||F<CR>

Table B.2-6. Recording patient information

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	8.1.1	Message type	P	1	Yes
2	8.1.2	Serial number	1	6	Yes
3	8.1.3	LIS number	101	5	No
6	8.1.6	Patient Name	ABC	30	No
9	8.1.9	Sex	M, F, U	1	No

Note.

- 1) Since only the first and second fields are mandatory, the patient information record is usually presented in the format "P|1 <CR>".
- 2) The total number of fields in the protocol is 35. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.
- 3) The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

B.2.4.3 Recording analysis information (O)**Note.**

The information in this section is provided in accordance with Section 9 of ASTM Protocol E1394-97. The information for each analysis should be presented using a record of this type.

Example:

O|1|1234567||^ALT^1:50|R<CR>

Table B.2-7. Recording analysis information

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	9.1.1	Message type	O	1	Yes

- 1) The first characters in the fifth field — "Analysis" — must be three characters "^", followed by the actual name of the analysis ("ALT" in the example above).
- 2) In the sixth field — "Priority" — the value "S" is used for the STAT sample, and the value "R" for the standard sample. The "R" value is mainly used.
- 3) The total number of fields in the protocol is 31. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.
- 4) The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

Field		E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1		10.1.1	Message type	R	1	Yes
2		10.1.2	Serial number	1	5	Yes
3	3.4	10.1.3	The studied indicator	ALT	10	Yes
	3.6		Dilution coefficient	1:50	10	No
4		10.1.4	Result	115.3	12	No

5	10.1.5	Unit measured	pg/mL	10	No
6	10.1.6	Reference range	0 to 200	30	No
7	10.1.7	Result label	L, H, N	1	No
13	10.1.13	Time to complete the analysis	YYYYMMDD HHMMSS	14	No

Note.

- 1) The first characters in the third field — "Analysis" — must be three characters "^", followed by the actual name of the analysis ("ALT" in the example above).
- 2) In the seventh field "Result label", the symbol "<" indicates a value below the linearity range, ">" indicates a value above this range. "L" means that the value is below the normal range, "H" means that the value is above the normal range, "N" indicates the normal value.
- 3) The total number of fields in the protocol is 14. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.
- 4) The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

**B.2.4.5 Recording the query for information (Q)****Note.**

The information in this section is provided in accordance with Section 12 of ASTM Protocol E1394-97. This record type is used to query the analysis information associated with a sample from the LIS system.

Example:

Q1^1234567||ALL|||||O<CR>

Table B.2-9. Recording the query for information

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	12.1.1	Message type	Q	1	Yes
2	12.1.2	Serial number	1	6	Yes

3	12.1.3	Sample ID	^1234567	22	Yes
5	12.1.5	All indicators under study	All	10	Yes
13	10.1.13	Status of information query	O	1	Yes

Note.

- 1) The first character in the third field, "Sample ID", must necessarily be a "^" character followed by the actual sample ID ("1234567" in the example above).
- 2) The total number of fields in the protocol is 13. The table above does not include all fields required by this software. The remaining fields can be added directly during data transfer; they will be automatically recognized by the software.
- 3) The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

B.2.4.6 Recording of the comments (C)

Note.

The information in this section is provided in accordance with Section 11 of ASTM Protocol E1394-97. It is intended to describe the results of the analysis.

Example:

C|1||1\2\3|I<CR>

Table B.2-10. Recording of the comments

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarity
1	11.1.1	Message type	C	1	Yes
2	11.2.1	Serial number	1	6	Yes
3	11.3.1	Source of the comment	I	1	Yes
4	11.4.1	Text of the comment	1\2\3	90	Yes
5	11.5.1	Type of the comment	I	1	Yes

B.2.4.7 Recording-end of the message (L)**Note.**

The information in this section is provided in accordance with Section 13 of ASTM Protocol E1394-97. Used in the last record of the transfer, indicates that the transfer has been completed.

Example:

L|1|N<CR>

Table B.2-11. Finalizing the message

Field	E1394	Field name according to ASTM protocol	Content	Maximum length	Necessarily
1	13.1.1	Message type	L	1	Yes
2	13.1.2	Serial number	1	6	Yes
3	13.1.3	Final code	N	1	Yes

**Note.**

The use of the <CR> delimiter at the end of message entries is mandatory. This delimiter should be added at the end of each message record to indicate the end of each record.

B.2.5 Example**B.2.5.1 Query of analysis details**

After placing the rack in the sample area and scanning the sample barcode, the software requests the corresponding sample analysis information from the LIS system by sending the message below.

Content of the message

```

--><ENQ>
<--<ACK>
--><STX>
<--<ACK>
-->H|^&||PSWD|MAGLUMI X8 ||||Lis||P|E1394-97|20100323<CR>
Q|I|^1234567||ALL|||||O<CR>
L|I|N<CR>
<--<ACK>
--><ETX>
<--<ACK>
--><EOT>
<--<ACK>

```

Table B.2-12. The value of the symbols

Code	Description	Corresponding ASCII code
-->	Sending by software	
<--	Getting the software	
<ENQ>	Query	0x05
<ACK>	Confirmation of response	0x06
<STX>	Start of the document	0x02
<ETX>	End of the document	0x03
<CR>	Start of the line	0x0D
<EOT>	Completing data transfer	0x04

Where:

```

H|^&||PSWD|MAGLUMI X8 ||||Lis||P|E1394-97|20100323<CR>
Q|I|^1234567||ALL|||||O<CR>
L|I|N<CR>

```

The above string is a message in ASTM E1394-97 format querying the reagent for the sample. The above text is a query for the reagent corresponding to the pattern 1234567.

Note.

- 1) In the above message, the <ACK> character is an acknowledgement returned by the LIS system. When this message is received, the appropriate <ACK> command must be sent in the above position. Otherwise, no connection to the LIS system will be registered in the software.
- 2) After the LIS system receives the message, it sends a reply message

with information about the analysis.

B.2.5.2 Response message with information about the analysis

After the LIS system receives the query from the software to send analysis data, a response message containing the relevant information is sent by the LIS system to the software.

Content of the message

```
<--<ENQ>
--><ACK>
<--<STX>
--><ACK>
<--H|^&||PSWD| MAGLUMI X8 ||||Lis||P|E1394-97|20180817<CR>
P|1<CR>
O|1|1234567||^CA125|R<CR>
O|2|1234567||^CA153|R<CR>
O|3|1234567||^CYFRA211|R<CR>
O|4|1234567||^FT3|R<CR>
O|5|1234567||^FT4|R<CR>
O|6|1234567||^T3|R<CR>
O|7|1234567||^TG|R<CR>
O|8|1234567||^TGA|R<CR>
L|1|N<CR>
--><ACK>
<--<ETX>
--><ACK>
<--<EOT>
--><ACK>
```

Where:

```

H|^&||PSWD| MAGLUMI X8 ||||Lis|P|E1394-97|20180817<CR>
P|1<CR>
O|1|1234567||^CA125|R<CR>
O|2|1234567||^CA153|R<CR>
O|3|1234567||^CYFRA211|R<CR>
O|4|1234567||^FT3|R<CR>
O|5|1234567||^FT4|R<CR>
O|6|1234567||^T3|R<CR>
O|7|1234567||^TG|R<CR>
O|8|1234567||^TGA|R<CR>
L|1|N<CR>

```

The above line is a response message in ASTM E1394 format sent by the LIS system and containing information about the analysis of the corresponding sample. The above text is the tests (CA125, CA153, CYFRA211, FT3, FT4, T3, TG and TGA) sent in response by the LIS system for sample 1234567.

Note.

In the above message, the <ACK> character is an acknowledgement returned by the LIS system. When this message is received, the appropriate <ACK> command must be sent in the above position. Otherwise, no connection to the LIS system will be registered in the software.

B.2.5.3 Sending test results

Content of the message

```

--><ENQ>
<--<ACK>
--><STX>
<--<ACK>
-> H|^&||PSWD| MAGLUMI X8 ||||Lis|P|E1394-97|20180817<CR>
P|1<CR>
O|1|1234567||^CYFRA211<CR>
R|1|^CYFRA211|0.8|ng/mL|0 to 7|N|||||20100326172956<CR>
L|1|N<CR>
<--<ACK>
--><ETX>
<--<ACK>
--><EOT>
<--<ACK>

```

Where:

```
H|\^&||PSWD| MAGLUMI X8 ||||Lis||P|E1394-97|20180817<CR>
P|I<CR>
O|I|1234567||^CYFRA211<CR>
R|I|^CYFRA211|0.8|ng/mL|0 to 7|N|||||20100326172956<CR>
L|I|N<CR>
```

The above line is an ASTM E1394-97 message with the results of the analysis sent to the LIS system. The above text is the CYFRA211 analysis results message for sample 1234567 that the software sends to the LIS system.

Note.

In the above message, the <ACK> character is an acknowledgement returned by the LIS system. When this message is received, the appropriate <ACK> command must be sent in the above position. Otherwise, no connection to the LIS system will be registered in the software.

B.3 HL7 Protocol

B.3.1 Description of the standard data transfer format

Format

<SB>CONTENT<EB><CR>.

Example

<SB>TEST<EB><CR>

Table. 8.1-1. Description of the standard data transfer

Code	Description	Corresponding ASCII code
<SB>	The initial character	0x0B
T	The letter T	0x54
E	The letter E	0x45
S	The letter S	0x53
T	The letter T	0x54
<EB>	End symbol	0x1C
<CR>	Carriage return symbol	0x0D

B.3.2 Description of the separators

Table B.3-2. Description of the separators

Separators	Description	Corresponding ASCII code
	Field separator	0x7C
—	Repeat separator	0x7E
^	Component separator	0x5E
&	The sub-component separator	0x26
\	Shielding symbol	0x5C

B.3.3 Description of the message types

Table B.3 -3. Types of the messages

Message label	Description
OML^O33	Sending analysis information in response
TSREQ	Query of analysis details
OUL^R22	Sending the analysis

	results
RSREQ	Query of analysis results

B.3.3.1 Response message with information about the analysis (OML^O33)

Note.

This message is used if the LIS system initiates sending orders in the analyzer or if the LIS system sends a response to the analyzer's orders.

Table B.3-4. The structure of the OML^O33 message

Message Segment	Value	Comment
MSH	The title of the message	
[PID]	Identification of the patient	Optional
SPM	Sample	
ORC	General order segment	
{OBR}	Research order	The recurring element

Note.

The total number of message segments in the protocol is 7. The list contains only the message segments used in this component. Other message segments can be added directly during data transfer; they will be automatically recognized by the system.

B.3.3.2 Analysis data query (TSREQ)

Note.

This message is used by the analyzer when requesting analysis data from the LIS system.

Table B.3-5. The structure of the TSREQ message

Message Segment	Value	Comment
MSH	The title of the message	
QPD	Defining query parameters	
RCP	Response Management Options	

B.3.3.3 Transmission results (OUL^R22)

Note.

This message is used if the analyzer initiates sending analysis results to the LIS system or sends a response message with the results of the query to the LIS system.

Table B.3-6. The structure of the OUL^R22 message

Message Segment	Value	Comment
MSH	The title of the message	
SPM	Sample	
ORC	General order segment	
{OBR}	Research order	The recurring element
OBX	Research/result	The recurring element
[NTE]}	Comments	Optional, a repeating element

**Note.**

The total number of message segments in the protocol is 11. The list contains only the message segments used in this component. Other message segments can be added directly during data transfer; they will be automatically recognized by the system.

B.3.3.4 Query Result (RSREQ)**Note.**

This message is used if the LIS system requests the analysis results from the analyzer.

Table B.3-7. RSREQ message

Message Segment	Value	Comment
MSH	The title of the message	
QPD	Defining query parameters	
RCP	Response Management Options	

B.3.4 Description of message fields**Note.**

The total number of message segments in the protocol is 13. The list contains only the message segments used in this system. Other message segments can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.1 Title (MSH)**Format:**

MSH|^~\&|Sender||Recipient||The time of sending the message||Message type|Confirmation code|P|2.5||NE|NE||UTF-8

Example:

MSH|^~\&|MAGLUMI||LIS||20151025105016||TSREQ|013329845|P|2.5||NE|NE||UTF-8

Table B.3-8. Description of the MSH message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	Field separator		1	Yes	Fixed value
2	Character encoding	^~&	4	Yes	Fixed value
3	Sender application	MAGLUMI	180	Yes	
5	Recipient application	LIS	180	Yes	
7	The time of sending the message	20151025105016	14	Yes	YYYYMMDDHHMMS S
9	Message type	TSREQ	7	Yes	OML^O33, OUL^R22, RSREQ or TSREQ
10	Control element ID/Message confirmation code	013329845	20	Yes	Unique ID
11	Processing ID	P	3	Yes	Fixed value
12	ID version	2.5	8	Yes	Fixed value
15	Confirmation of receipt of the message	NE	2	Yes	Fixed value
16	Confirmation from the application	NE	2	Yes	Fixed value
18	Symbol set	UTF-8	16	Yes	Fixed value

**Note.**

The total number of fields in the protocol is 18. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.2 Patient Identification (PID)**Format:**

PID|1||Patient's ID||Patient's full name||yyyyMMdd(Date of birth)|Gender
(Enumeration)|||||||||Age^The unit of measurement of age

Example:

PID|1||10248734||Zhang Jianguo||20110101|1|||||||8^1

Table B.3-9. Description of the PID message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	№	1	1	Yes	Fixed value
3	Patient's ID	10248734	20		
5	Patient Name	Zhang Jianguo	50	Yes	
7	Date or time of birth	20110101	8		YYYYMMDD
8	Sex	1	1	Yes	M: male, F: female. U: unknown
38	38.1	Age	8	Yes	Positive integer
	38.2	Unit of measure of age	1	Yes	1: year, 2: month, 3: day, 4: hour

Note.

The total number of fields in the protocol is 38. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.3 Sample (SPM)

Format:

SPM||ID of sample ^ ^ID of control material^ID of sample rack^Sample rack position|| Sample type|||||Sample role|||||

Example:

(Sample) SPM|1|034581229010^101^^A001^2||0|||||P|||||

(QC) SPM|1|QC material name^QC lot number^#2020030301#^C01^4||0|||||Q|||||

Table B.3-10. Description of the SPM message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	№	1	1	Yes	Fixed value
2	2.1	Sample ID	034581229010	22	Yes
	2.2	sample no. (LIS number)	101	5	
	2.3	ID of the control material	#2020030301#	22	Yes
	2.4	ID of the sample rack	A001/C01	22	

2.5	Position of the sample rack	2/4	5		
4	Sample type	1/6	1	Yes	0: serum 1: blood plasma 2: urine 3: hydrothorax and ascites 4: cerebrospinal fluid 5: Other
11	Role of the sample	P/Q	1	Yes	P: Patient sample Q: Control materials

**Note.**

1. The total number of fields in the protocol is 18. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.
2. When the sample is used as QC material, Field 2 is the name of the QC material and Field 2.2 is the QC lot number.

B.3.4.4 Order total segment (ORC)**Format:**

ORC|Status of the order|||||Priority|||Time of the order|||

Example:

ORC|NW|||||0||201909131511433|||

Table B.3-11. Description of the ORC message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	Status of the order	NW	2	Yes	NW: new physician's order RF: repeat order query
7	Priority	0	1	Yes	0: General 1: Quick check
9	Time of the order	201909131511433	14	Yes	YYMMDDHHMMS S

**Note.**

The total number of fields in the protocol is 12. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.5 Research order (OBR)

Format:

OBR||||| Code of Analysis^Exploration Coefficient.

Example: OBR|1||||FT3^1:10

Table B.3-12. Description of the OBR message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	№	1	4	Yes	Starts at 1
4	4.1	Anal. code	FT3	Yes	Reagent channel no.
	4.2	Dilution coefficient	1:10		

Note.

1. The total number of fields in the protocol is 4. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.
2. Specify whether to send a dilution factor on the [Network Connection] screen.

B.3.4.6 Result (OBX)

Format:

OBX|||||Analysis code||Analysis result||Units of measurement||Reference range||Anomaly label|||F||Begin time of analysis|| |Module name|Equipment name|Analysis end time||Reference label|Message exchange type

Example:

OBX|1||FT4||20.9|mmol/L|8.9-17.2|H||F||20191025104016|||Module1^Maglumi|20191025104816|0|0

Table B.3-13. Description of the OBX message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	№	1	4	Yes	Starts at 1
3	Anal. code	FT4	10	Yes	Abbreviated name of the reagent
5.1	Quantitative result	20.9	10	Yes	

6	Unit of measure	2	10		Units of measure for reagent analysis results
7	Reference range	8,9–17,2	100		Reference range of values for reagent assay results
8	Abnormal value label	H	6		H: above normal HH: above the linearity range L: below normal LL: below the linearity range N: no Label
11	Result Status	F	1	Yes	Fixed value, indicates that the result is final
14	Analysis start time	20191025104016	14	Yes	YYYYMMDDHHMMS S
18	18.1 Module Name	Module1	20	Yes	
	18.2 Equipment name	Maglumi	20	Yes	
19	Analysis end time	20191025104816	14	Yes	YYYYMMDDHHMMS S
20	Retest label	0	1	Yes	0: Result of standard analysis 1: Result of retest
21	Type of messaging	0	1	Yes	0: Automatic messaging 1: Manual messaging

**Note.**

The total number of fields in the protocol is 21. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.7 Notes and Comments (NTE)

Format:

NTE|№

Example:

NTE|1||1~2

Table B.3-14. Description of the NTE message segment

Position	Field	Content	Maximum length	Necessarily	Comment
----------	-------	---------	----------------	-------------	---------

1	№	1	4	Yes	Starts at 1
3.1	Code	1–2	10		

**Note.**

The total number of fields in the protocol is 3. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.8 Query Parameters (QPD)**Format:**

QPD|Query Name||Sample ID^Sample Rack ID^Sample Rack Position|Priority

Example:

QPD|RSREQ||123456^NA001^2|0

Table B.3-15. Description of the QPD message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	Name of the query message	RSREQ	7	Yes	RSREQ or TSREQ
3	3.1	Sample ID	123456	Yes	
	3.3	ID of the sample rack	NA001		
	3.4	Position of the sample rack	2		
4	Priority	0		Yes	0: General 1: STAT

**Note.**

The total number of fields in the protocol is 4. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.4.9 Response control (RCP)**Format:**

RCP|Query priority||| Response form

Example:

RCP||R|

Table B.3-16. Description of the RCP message segment

Position	Field	Content	Maximum length	Necessarily	Comment
1	Query priority	I	4	Yes	Fixed value, indicates the urgency of the query
3	Response Form	R	10	Yes	Fixed value, indicates that responses are given in real time

Note.

The total number of fields in the protocol is 3. The list contains only the message segments used in this component. Other fields can be added directly during data transfer; they will be automatically recognized by the system.

B.3.5 Example

B.3.5.1 Analysis data query (

Actual message:

```
--> MSH|^~\&|MAGLUMI X8||Lis||20200311140350||TSREQ|2|P|2.5|||NE|NE||UTF-8
--> QPD|TSREQ||1234567^^A006^2|0
--> RCP|I||R|
```

Description of the messaging process

Code	Description
-->	Sending by software

The analyzer sends the above message to the LIS system using the HL7 data transfer protocol to request information about the analysis of the corresponding sample. The above message text performs the query for sample analysis information 1234567. Details on the meaning of the fields are provided in Appendix B, Section 3.4.

B.3.5.2 Response message with information about the analysis

Actual message:

```
<-- MSH|^~\&|Lis||MAGLUMI X8||20200311140350||OML^O33|2|P|2.5|||NE|NE||UTF-8
<-- PID|1||34||zhangsan||19960101|Male|||||||||||||||||^ years old
```

```

← SPM|1|1234567^^|P
← ORC|NW|||N
← OBR|1||FT3^
← OBR|2||FT4^

```

Description of the messaging process:

Code	Description
←	Getting the software

The LIS system sends the above message using the HL7 data transfer protocol to send the analysis information of the corresponding sample in response. The text of the above message executes the sending of FT3 and FT4 analysis details for sample 1234567 in response. This message is the opposite of the message in Section 3.5.1 of Section B, which describes the process for querying analysis details. Details on the meaning of the fields are provided in Appendix B, Section 3.4.

B.3.5.3 Result of messaging**Actual message:**

```

→ MSH|^~&|MAGLUMI X8|Lis|20200311140933|OUL^R22|4|P|2.5||NE|NE|UTF-8
→ SPM|1|1234567^^^A006^2|0||P|
→ ORC|NW|||0|20200311140350||
→ OBR|1||FT4
→ OBX|1||FT4||9.2|pg/mL|8.90-45.00|HH||F||20200311140702||Module1^MAGLUMI|20200311140933|
0|0
→ NTE|1

```

Description of the messaging process:

Code	Description
-->	Sending by software

The analyzer sends the above message using the HL7 data transfer protocol to send the analysis results to the LIS system. The above text is the FT4 analysis results message for sample 1234567 that the analyzer sends to the LIS system. Details on the meaning of the fields are given in chapter B.3.4.

B.3.5.4 Query Result**Actual message:**

```

<-- MSH|^~^&|Lis|MAGLUMI X8||20200311142249||RSREQ|5|P|2.5||NE|NE||UTF-8
<-- QPD|RSREQ||1234567^|N
<-- RCP|||R|

```

Description of the messaging process:

Code	Description
<--	Getting the software

The LIS system sends the above message to the analyzer using the HL7 data transfer protocol to query the analysis result. The text of the above message is querying the result of sample 1234567. Details on the meaning of the fields are given in chapter B.3.4.

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Snibe, **MAGLUMI**, **Blossays**, **Biolumi**, **Preaccu**, **SnibeLinker**, **Molecision** and **SnibeLis** are registered trademarks or trademarks of Shenzhen New Industries Biomedical Engineering Co., Ltd.