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CERTIFICATE

中国国际贸易促进委员会暨中国国际商会
China Council for the Promotion of International Trade is China Chamber of International Commerce

01701844

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of SIEMENS SHENZHEN MAGNETIC RESONANCE LTD. on the annexed DOCUMENT is genuine.

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证书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Document for submission


Hereby certify that the following document has been developed, prepared, and released by the manufacturer of the medical device.

User Manual supplement for Russia – Angiography system ARTIS one with accessories

Sincerely yours,
Siemens Shenzhen Magnetic Resonance Ltd.

NAME: Yin Li Nuan

POSITION: Regulatory Affairs & Standard Engineer

SIGN  2023.10.26

STAMP

Siemens Shenzhen Magnetic Resonance Ltd.

No.32 Gaoxin C. Ave
2nd, Hi-Tech Industrial Park
Shenzhen 518057,
PEOPLE'S REPUBLIC OF CHINA

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User Manual addendum for Russia

Angiography system ARTIS one with accessories

1. The name of the Medical Device

Angiography system ARTIS one with accessories.

*The text of the document may also contain the following names for the above-mentioned medical device:
"ARTIS one system", "ARTIS one".*

2. Manufacturing information of the Medical Device

Legal manufacturer:

Siemens Shenzhen Magnetic Resonance Ltd.,
Siemens MRI Center,
Gaoxin C. Ave., 2nd Hi-Tech Industrial Park,
518057 Shenzhen, China

Manufacturing sites:

Siemens Shenzhen Magnetic Resonance Ltd.,
Siemens MRI Center,
Gaoxin C. Ave., 2nd Hi-Tech Industrial Park,
518057 Shenzhen, China

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubininskaja street 96,
115093 Moscow, Russia
Tel.: 8 800 550 50 22
www.siemens-healthineers.com/ru/

3. Designation of the medical device

3.1. Intended use

ARTIS one, as a family of dedicated angiography systems generating fluoroscopic and radiographic X-ray images, is intended for diagnostic imaging, interventional and minimally invasive therapy.

Its use is indicated for image guidance in non-vascular, vascular, cardio-vascular, neuro-vascular, minimal-surgical, and surgical procedures. The ARTIS one system includes also software options which allow the reconstruction of two-dimensional images acquired with rotational angiography into a three-dimensional image format, as well as options supporting diagnosis, planning, guidance and treatment follow-up.

3.2. Indications

In general, all clinical indications of angiography- and fluoroscopic-based procedures within the intended use are applicable for this device.

ARTIS one is an angiography system developed for diagnostic imaging and interventional procedures including, but not limited to, pediatric and obese patients.

Procedures that can be performed with the ARTIS one include cardiac angiography, neuro-angiography, general angiography, rotational angiography, multipurpose angiography and whole body radiographic/fluoroscopic procedures as well as procedures next to the table for i.e. patient extremities. Additional procedures that can be performed include angiography in the operating room, image guided surgery by X-ray, by image fusion, and by navigation systems. The examination table as an integrated part of the system can be used for X-ray imaging, surgery and interventions. The ARTIS one can also support the acquisition of position triggered imaging for spatial data synthesis.

3.3. Contraindications

For this product there are currently no known contra-indications. However, contraindications for angiography- and fluoroscopy- based procedures apply.

The final decision for use of the medical device in a certain application is made by the physician based on his/her medical knowledge and the risks of the cases involved.

3.4. Undesirable side effects

Deterministic and stochastic effects of exposure to ionizing radiation (X-rays) are part of the physician's medical education. The necessity of an X-ray examination in relation to these risks must be considered by the physician when prescribing or performing such a type of examination.

For radiation emitting products, known deterministic effects of exposure to X-rays are erythema, cataracts, permanent epilation and delayed skin necrosis.

For radiation emitting products, known stochastic effects of exposure to X-rays are increased risks of developing cancer and hereditary diseases.

These stochastic dose effects have an increased probability of occurrence with increased dose. The related risk can be minimized by means of keeping the X-ray exposure to the necessary minimum and by means of careful and systematic application of protection measures for the operator and the clinical staff.

3.5. Patient target group(s)

The system can be used on all patients, from newborn to geriatric. Patient weight is limited to the specification of the patient table.

3.6. Intended users

- Operator profile: The usage of the system described in the Operator Manual requires specific technical and medical knowledge and skills regarding, at a minimum, radiation protection, safety procedures and patient safety.

People using, moving, working with the system must have acquired such knowledge and skills during their curriculum.

- Equipment training: Application training is delivered with the equipment according to the handover contract. It is mandatory to follow such application training delivered by Siemens Healthineers Representative before any use of the system.

The follow-up training, which is necessary due to change of personnel, is in the responsibility of the operator of the system. Any additional training can be requested from Siemens Healthineers.

3.7. Environmental conditions for operation

General Room Specifications for: Examination room, Control room, Equipment room

Temperature range: 15 – 30°C, 22°C is recommended;

Relative humidity: 20 – 75% without condensation;

Max. temperature gradient: 5 K/h;

Barometric pressure: 70 – 106 kPa.

3.8. Environmental conditions for transport and storage

Temperature range: -20 – 70°C;


Relative humidity: 10 – 95% without condensation;

Barometric pressure: 70 – 106 kPa.

4. System's overview

4.1. ARTIS one system overview

The individual components of your system vary depending on the configuration. In the following you see an overview of a typical configuration.

Example


- (1) Floor stand with C-arm, X-ray tube assembly and FD
- (2) Display ceiling suspension (DCS) with 30" Display, optionally with additional LC display
- (3) Patient table
- (4) Control console for controlling the stand, patient table, and imaging system
- (5) Handswitch
- (6) Footswitch for releasing radiation

Note: The cover of the system includes an official brand name of the system «ARTIS one». Also the Cover of the system includes the following information:

➤ Edition X

This additional information isn't a system's name and this information isn't indicated on the main system's labeling, User documentation etc. This additional information isn't part of the system's name and intends to be informative of the technology features of the system.

4.2. The floor stand


- (1) C-arm
- (2) FD slide
- (3) Flat detector (FD)
- (4) X-ray tube assembly with primary collimator
- (5) Swivel base
- (6) Unlocking pushbutton
- (7) Stand column

The floor stand is mounted on a swivel base. The floor stand can be moved between working and patient transfer positions on the swivel base.

The floor stand can be moved into various positions relative to the patient. Oblique projections are possible by rotational and orbital movement (angulations) of the C-arm. The X-ray tube assembly with primary collimator and the flat detector (FD) are attached to the C-arm. The FD is mounted on a lifting device for a variable SID.

4.3. Patient table

- (1) Tabletop with mattress
- (2) Accessory rails for control modules and other accessories
- (3) Sockets for injector and country specific power outlet
- (4) Equipotentiality: Connectors for equipotential equalization
- (5) Pedal for releasing the table rotation lock

Connector for the footswitch

The ARTIS one system is equipped with a floor-mounted patient table.

The table height and tabletop longitudinal and transverse positions are continuously adjustable. The patient table can also be rotated.


Rails for control modules and other accessories are attached to the table. Devices with equipotential equalization can be connected to the table (4).

4.4. Control console

Your system is equipped with one or two control consoles.

The control console consists of one grouped or several individual control modules.

There are four basic types of control modules.

Control console (example)


- (1) Table (Table Control Module TCM)
- (2) Stand/C-arm (Stand Control Module SCM)
- (3) Collimation and filter (Collimator Control Module CCM)
- (4) Imaging system (Table-Side Control TSC)
- (5) Handle for attaching the control module

The control modules are equipped with a joystick and several keys.

4.5. Footswitch

Footswitches with four pedals are connected to the system for releasing fluoroscopy, acquisition, or other configurable functions.

Additionally, a wireless footswitch may be available in the examination room.

Example configuration

- (1) Fluoroscopy
- (2) Acquisition
- (3) Alternative acquisition*
- (4) Tabletop brake*
- (5) Safety bracket

* The function of the footswitch pedals can be configured by Customer Service.

4.6. Wireless footswitch

The wireless footswitch has the same functions as the standard footswitch without being connected by cable. The functions are activated via an RF transmitter.


The immanent advantages are easier cleaning and no tripping over cables.

The wireless footswitch is active in a range of approximately 5 m (16 ft) around the table/stand, where the receiver is mounted. If the footswitch is in a wider range of the table, it may lose contact with the receiver and will then turn into an inactive mode.

4.7. Display Ceiling Suspension (DCS)

Depending on system configuration, the DCS carries one or two displays. Displays of sizes 30" or 55" can be installed.

Example

Using the handle, the DCS can be moved, rotated, and adjusted in height. An indicator light is integrated on the DCS boom above the display.

4.8. Mode of application:

Before starting an examination, get used to:

- Switch-on the system;
- Perform safety checks;
- Prepare the stand and table by performing unit movements;
- Adjust the image format, collimation, and filters;
- Insert or remove grid;
- View and configure the screen layout;
- Operate the imaging system in the examination room.

5. Classification of the Medical Device

5.1. Type of the medical device in accordance with the nomenclature for classification

Classification according to Annex IX, MDD

Class

IIb

Classification in accordance with Russia Nomenclature classification:

191140

5.2. Class of potential Risk of the medical device application

The Medical Device is intended for transmission of energy in the form of ionizing radiation, it belong to the Class 2b.

5.3. Type on contact with the human body

The Medical Device is non-sterile, nonexpendable.

5.4. Software classification

The Medical Device's Software safety Class according to IEC 62304 is Class B.

SW version: VA11; date of release: 03/2022.

5.5. Applied parts

According to IEC 60601-1 Cl. 7.2.10 and 8.3 the following parts are applied parts of **type B**:

- Tabletop (upper surface);
- Patient mattress;
- Table accessories, intended to support patient (e.g. extensions, head holder etc.):
 - Head Wedge;
 - Head position (prone - adult);
 - Head pad, deep;
 - Head pad, flat;
 - Handle;
 - Hand grips with supports;
 - Arm rest, cpl;
 - Arm rest;
 - Handgrip;
 - Arm holder, left and right, 1 pair;
 - Set of body straps.

The system is not classified as defibrillator-proof applied part (IEC 60601-1 Clause 8.5.5).

Note: following components are identified as parts that can come into contact with PATIENT but not defined as APPLIED PARTS, subjected to the requirements for APPLIED PARTS, except for Clause 7.2.10: detector cover, X-ray tube incl. housing, collimator cover, C-arm (inner surfaces), Control modules.

5.6. Class of the medical device according to IEC 60601-1

According to clause 6.2 of IEC 60601-1:2005, the Medical Device is the Class 1.

5.7. Operating mode

Operating conditions according to IEC 60601-1.

System: Continuous mode of operation

Patient table: 6 min/hour

5.8. Protection against water ingress is classified per IEC 60529

IPX2: Patient table;

IPX4: Control modules, handswitch;

IPX8: Footswitches;

IPX0: Rest system, including: generator, flat panel detector, control console, DAP chamber, UPS and Raysafe Real-time Radiation Dosimeter.

5.9. AP/APG

The system is not intended for use with flammable anaesthetics (Category AP) or flammable anaesthetics with oxidants (Category APG) and in oxygen rich environment.

6. Principles of the Medical Device operation

The ARTIS one system consists of the acquisition system, the imaging system, and supporting devices like the generator.

The components of the system are installed in the examination room and in the control room. Supporting devices may be installed in different rooms.

ARTIS one is a fully digital angiography system equipped with flat detector technology.

Depending on the system configuration, an ARTIS one is suitable for innovative cardio-vascular examinations and interventions, universal angiography procedures and interventional procedures.

7. Technical Data: Specifications and Characteristics


7.1. Technical Data

Value statements All mechanical values listed below have a tolerance of $\pm 5\%$.

Measuring functions $\pm (5\% + 1 \text{ last digit})$

Accuracy of radiation output $\pm 5\%$

System view ARTIS one in mm

Room height for Large display RH 2900 - 3100 (114.17 - 122.05")


ARTIS one stand positions

Transfer position

Head-end position (e.g. for PCI)

Left-side position


Right-side, table-rotated position
(e.g. for pacemaker implantations)

Peri-position

Left-side, OR position (e.g. for AAA treatment)

Extended parking position*

OR position, right-side*

Additional transfer position*

7.1.1. ARTIS one

C-arm positioning	The C-arm positioning (RAO, LAO, cran, caud) depends on the installation, table top, table position and patient size. Please refer to the values displayed on the data display.
Variable focal spot-to-detector distance	90 cm - 120 cm
Focus-to-isocenter distance	75 cm
Isocenter-to-floor distance	107 cm
Patient coverage (free floating tabletop, minimum without repositioning)	210 cm with table movement
C-arm depth	92.5 cm
Stand rotation	motorized programmable positioning
Angulation speed	variable rotation up to 25°/s with LAO/RAO and 25°/s with CRAN/CAUD; variable rotation, 3D up to 60°/s
Longitudinal C-arm movement	motorized up to 250 mm/s
Transversal C-arm movement	motorized up to 150 mm/s
C-arm oblique projections ¹⁾	± 130° LAO/RAO and + 55°/- 45° CRAN/CAUD at 0° head-end C-arm position; + 81° LAO to - 59° RAO and + 48° CRAN to - 53° CAUD at 35° left-side C-arm position
Weight ²⁾ :	Stand: 570 kg

¹⁾Maximum angulations depend on stand position, table position and patient size.

²⁾ All weight specifications are for reference only. Depending on the system configuration, the resulting floor loads will be different.

7.1.2. Patient table

Table height	75 cm - 110 cm
Table length	284.5 cm
Table rotation	±120°
Weight ¹⁾	255 kg
Maximum table load	• 250 kg patient weight on tabletop • 100 kg table accessories on rails – 60 kg on side rails – 40 kg on foot rail • 60 kg cardiopulmonary resuscitation (CPR)
Attenuation equivalent according to 21CFR1020.32 / EN60601-1-3 (100 kV @ 3.6 mm HVL)	Tabletop with mattress: ≤ 2.0 mm Al

Manual longitudinal travel	125 cm
Manual transverse travel	±17,5 cm
Duty cycle max. 10%	The lifting motor is specified for usage of max. 60 s/10min.
Lift speed	5 cm/s

¹⁾ All weight specifications are for reference only. Depending on the system configuration, the resulting floor loads will be different.

7.1.3. Radiation generation

X-ray generator

Model: Polydoros A100 Plus

Manufacturer: Siemens Healthcare GmbH

Mode of operation	Continuous operation with intermittent, short time loading (duty cycle controlled by heat unit calculator)
Nominal x-ray tube voltage @ highest x-ray tube current	100 kV @ 1000 mA
Highest x-ray tube current @ Maximum x-ray tube voltage	800 mA @125 kV
Highest electrical output power	100 kW, 100 kV @1000 mA
Nominal electrical output power	100 kW, 1000 mA @100 kV, 100 ms
Lowest current-time product	48 µAs (15 mA @3.2 ms) @ 40 - 125 kV (in combination with ARTIS one)
X-ray tube voltage	40 - 125 kV
X-ray tube current	0.5 - 1000 mA
Pulse time	2 ms - 800 ms
Pulse frequency	0.5 - 100 Hz
Accuracy	X-ray tube voltage (kV) : ±5% X-ray tube current (mA) : ±10% pulse time (ms): ±(10% +1 ms)

System including MEGALIX® Cat Plus 125/40/90-GW X-ray tube assembly

Model: MEGALIX Cat Plus 125/40/90-GW

Manufacturer: Siemens Healthcare GmbH

Automatic exposure control	
X-ray tube voltage	40 - 125 kV
X-ray tube current	15 - 800 mA (depends on focal spot)
Pulse time	2 - 800 ms
Anode track speed	47 m/s
Copper filter	0 / 0.1 / 0.2 / 0.3 / 0.6 / 0.9 mm Cu
Reproducibility of the Radiation output	variation coefficient: < 0,05
Automatic Control system	repetition accuracy (dose output) ≤ 2%
Accuracy dose area product	± 35 %
Accuracy reference air kerma	± 35 % (IEC 60601-2-43)
Radiography (incl. Cineradiography)	(large focus F2)

Nominal x-ray tube voltage @ highest x-ray tube current	1000 mA @ 100 kV (100 kW)
Highest x-ray tube current @ Maximum x-ray tube voltage	125 kV @ 800 mA (100 kW)
Highest electrical output power	100 kW (100 kV @ 1000 mA)
Nominal electrical output power	100 kW (100 kV @ 1000 mA)
Lowest current-time product @ x-ray tube voltage range	48 μ As (15 mA @ 3.2 ms) @ 40 – 125 kV
Radioscopy (Fluoroscopy)	(small focus F1)
Nominal x-ray tube voltage @ highest x-ray tube current	125 kV @ 250 mA
Highest x-ray tube current @ Maximum x-ray tube voltage	250 mA @ 125 kV
Highest electrical output power	31,25 kW (125 kV @ 250 mA)
Lowest current-time product @ x-ray tube voltage range	48 μ As (15 mA @ 3.2 ms) @ 40 - 125 kV
Accuracy	x-ray tube voltage (kV) : $\pm 5\%$ x-ray tube current (mA) : $\pm 15\%$ pulse time (ms) : $\pm (10\% + 1 \text{ ms})$

X-ray tube assembly

Quality equivalent filtration $\geq 1.5 \text{ mm Al}$ at 75 kV/HVL 2.7 mm Al (at use of aluminum collimator cover)

Collimator

Max. field size (SID = 90 cm)	287 mm x 261 mm
Min. field size (SID = 100 cm)	1 cm x 1 cm
Quality equivalent filtration	0.9 mm Al at 70 kV/HVL 2.7 mm Al
CAREfilter	Additional filtration: 0 / 0.1 / 0.2 / 0.3 / 0.6 / 0.9 mm Cu Note: 0.2 mm Cu additional filtration is more than 1 mm Al + 0.1 mm Cu required for pediatric applications.

Dose measuring chamber for collimator

Model: KermaX plus 120-160 EAS

Manufacturer: IBA Dosimetry GmbH

Type	For measuring and displaying the dose area product in $\mu\text{Gy} \times \text{m}^2$ with measuring chamber in the primary collimator. Measuring unit and display.
Quality equivalent filtration	$\geq 0.2 \text{ mm Al}$ at 75 kV/HVL 2.7 mm Al

7.1.4. Image acquisition

Imaging size	287 mm x 261 mm
Image display matrix	1560 x 1420 pixels
X-ray conversion technology	a-Si with CsI scintillator
Digitization depth	16-bit
Pixel pitch	184 μ m
Nyquist frequency	2.7 lp/mm
DOE (detective quantum efficiency) (RQA5)	0.5 lp/mm at 2 μ Gy: typical 60% 1 lp/mm at 2 μ Gy: typical 52% 1.5 lp/mm at 2 μ Gy: typical 46% 2.7 lp/mm at 2 μ Gy: typical 21%
MTF (modulation transfer function) at 1 lp/mm	0.5 lp/mm: 80% typical (according to IEC 62220) 1 lp/mm: 59% typical (according to IEC 62220) 1.5 lp/mm: 41% typical (according to IEC 62220) 2.7 lp/mm: 17% typical (according to IEC 62220)
Signal to electronic noise ratio (SENR)	min. 8 dB, typical 11 dB, max. 13 dB (RQA5, 1x1 binning, high gain)
Grid	Removable; Pb 15/80; fo = 105 cm
Collision protection - attenuation equivalent	0.15 mm Al / HVL 3.6 mm Al, 100 kV
Entrance of detector	0.38 mm Al / HVL 3.6 mm Al, 100 kV

7.1.5. Image display

The performance required to present the images for diagnostic purpose is as follows:

IEC 60601-2-43:2019	Diagonal screen measurement: $\geq 21"$ Resolution: $\geq 1200 \times 1600$ Typical luminance: ≥ 420 cd/m ² Contrast ratio: $\geq 1500:1$
---------------------	--

7.1.6. Footswitches

Data for wireless footswitches:

Active range	5 m (16 ft) around the receiver mounted to the table/stand.
Charging time	6 h
Operating time with full batteries	60 h
Compliant Bluetooth Specification	Bluetooth specification Ver 2.0 + EDR
Hopping frequency	2402 MHz - 2480 MHz F = 2402 + k MHz, k=0..78
Modulation method	GFSK, 1 Mbps
Effective radiated power	0 dBm (typ.)

Wireless footswitch Charger:

Dimensions (L x W x H)	160 x 140 x 87 mm
------------------------	-------------------

Input current	100 V – 400 mA; 240V - 200 mA
Output current	12 V; 1500 mA
Power	18W

Data for all footswitches:

Quantity of pedal	4
Press force	25 N

7.1.7. Injectors

At present, the following injectors have been approved:

- ACIST CVI
- Illumena
- MEDRAD Avanta
- MEDRAD ARTERION
- Medtron_Accutron_HP_D
- Medtron_Accutron_HP

At present, the following systems have been approved for IVUS:

- Volcano IVUS

For operation, technical description, models, and technical data, please see the documentation supplied by the manufacturer.

7.1.8. Radialis armboard

Radialis armboard	Weight	2.5 kg
	Maximum load	10 kg
L shaped board	Weight	0,8 kg
	Maximum load	0.5 kg
Attenuation equivalent	Radialis armboard with cushion: ≤ 2.0 mm Al (according to 21CFR1020.30 / EN60601-1-3 (100 kV @ 3.6 mm HVL))	

7.1.9. Installation data

Line voltage connection, 3-phase current

X-ray generator

Rated voltage 3Ph $\pm 10\%$	380 V	400 V	420 V	440 V	460 V	480 V
Line resistance of generator at 100 kW	$\leq 0,08 \Omega$	$\leq 0,09 \Omega$	$\leq 0,09 \Omega$	$\leq 0,10 \Omega$	$\leq 0,11 \Omega$	$\leq 0,12 \Omega$
Line resistance of generator at 80 kW	$\leq 0,10 \Omega$	$\leq 0,11 \Omega$	$\leq 0,12 \Omega$	$\leq 0,14 \Omega$	$\leq 0,15 \Omega$	$\leq 0,16 \Omega$
Frequency ± 1 Hz	50/60 Hz					60 Hz
Fuse internal	50 A					
Fuse external	63 A slow-blow					
Power consumption	8.0 kVA for fluoro (long-time loading >5 min) 160 kVA for acquisition (momentary peak loading <5 sec)					

System control cabinet

Rated voltage 3Ph $\pm 10\%$	380 V	400 V	420 V	440 V	460 V	480 V
Frequency ± 1 Hz	50/60 Hz					60 Hz
Fuse internal	8 A					
Fuse external	25 A					
Power consumption	max. 7.2 kVA (long-time loading >5 min)					

To avoid the risk of electric shock, this equipment must only be connected to supply mains with protective earth.

7.1.10. Ambient conditions for operation

General Room Specifications for:

Examination room, Control room, Equipment room

Temperature range	+ 15 °C to + 30 °C (recommended temp. 22 °C)
Relative humidity	20 – 75 % without condensation
Max. temperature gradient	5 °C/h
Barometric pressure	70 – 106 kPa

Imaging system

For climatic conditions, see table General Room Specifications.

Air flow	156 m ³ /h
max. Noise level	53 dB (A)

Generator

For climatic conditions, see table General Room Specifications.

Air flow	160 m ³ /h
max. Noise level	55 dB (A)

System control cabinet

For climatic conditions, see table General Room Specifications.

Air flow	295 m ³ /h
max. Noise level	55 dB (A)

Cooling unit

Cooling air:	+ 15 °C to + 30 °C (frost-free room)
Air flow:	950 m ³ /h
Max. noise level	55 dB (A) at 50 Hz; 59 dB (A) at 60 Hz

Stand

Mechanical impact	max. 10 g/16 ms
Vibrations	max. 0.1 g / 10 – 200 Hz
Max. noise level	< 55 dB (A)

7.1.11. Ambient conditions for transport and storage

The system will be stored and installed by Siemens service. The system will be fixed installed and bolted in the respective room. Therefore, local transport of the user is not possible.

Temperature range	- 20 °C to + 70 °C
Relative humidity	10 – 95 % without condensation
Barometric pressure	70 – 106 kPa

7.1.12. Connection of ARTIS by NETWORK/DATA coupling to other equipment

If ARTIS one is connected to other equipment beyond the control of Siemens Healthineers by NETWORK/DATA coupling, the following has to be considered:

- The characteristics of the NETWORK/DATA coupling necessary for ARTIS one are described by the "DICOM Conformance Statement" document (available via Siemens Internet).
- Cautions and warnings about possible failures of the NETWORK/DATA coupling are

described at the appropriate chapters of the operator manual.

The responsible customers' organization has to consider that:

- Connection of ARTIS one to a NETWORK/DATA coupling that includes other equipment could result in previously unidentified risks to patients, operators or third parties.
- These risks should be identified, analyzed, evaluated and controlled by the responsible customers' organization.
- Subsequent changes to the NETWORK/DATA coupling could introduce new risks and require additional analysis. These changes to the NETWORK/DATA coupling include:
 - Changes in NETWORK/DATA coupling configuration
 - Connection of additional items to the NETWORK/DATA coupling
 - Disconnecting items from the NETWORK/DATA coupling
 - Update or upgrade of equipment connected to the NETWORK/DATA coupling


7.2. Scattered Radiation in the significant zone of occupancy

In the significant zone of occupancy (main operation area) diagnosis and therapy are performed.


Type of examination: Angiography using the femoral puncture technique

Position and extent of the significant zone of occupancy

Significant zone of occupancy: Patient support horizontal, X-ray tube below patient support

Angulation:
0° RAO
0° CRAN / CAUD

■ Significant zone of occupancy Dimensions in cm

7.2.1. Measuring conditions

Measuring conditions of scattered radiation in the significant zone of occupancy according to EN / IEC 60601-1-3

Angulation	0° RAO, 0° CRAN/CAUD
Continuous fluoroscopy	125 kV / 3.0 mA
Automatic control	Fluoro CARD
Total beam filtration	>2.5 mm Al + 0.2 mm Cu with continuous fluoroscopy (in BYPASS FLUORO mode only) >2.5 mm Al + 0.6 mm Cu with automatic control (0.2 mm Cu corresponds to 7.1 mm AL)
Field size on the FD collision protection	18 cm x 18 cm
Water phantom	25 cm x 25 cm x 25 cm lying on the patient table with mattress
Center of water phantom	on central X-ray beam axis

1. Measurement without radiation protection equipment



2. Measurement with mobile lower body radiation protection (lead equivalent value Pb=0.5 mm) and upper body radiation protection (lead equivalent value Pb=0.5 mm)

The lower body radiation protection is located at the patient table in front of the significant zone of occupancy.

7.2.2. Air kerma rates


Measurement	Conditions
A	with radiation protection, continuous fluoroscopy (in BYPASS FLUORO mode only)
B	with radiation protection, with Automatic Control Exam Set "FL Card": 70.0 kV, kV Dose 109 kV, 40 nGy/p, 5 ms, 0,2 Cu - 0,6 Cu, 15 pulse/s
C	without radiation protection, continuous fluoroscopy
D	without radiation protection, with Automatic Control Exam Set "FL Card": 70.0 kV, kV Dose 109 kV, 40 nGy/p, 5 ms, 0,2 Cu - 0,6 Cu, 15 pulse/s

Height above floor


7.3. Isokerma maps of stray radiation


100 cm, vertical position

100 cm, horizontal position

150 cm, horizontal position

7.4. Dosimetric information

Interventional reference point	60 cm
Radiation field size	See Table 1a and 1b See Table 3 for dose factors of other zoom formats.
Pulse frequency	See Table 1a and 1b See Table 2 for dose factors of other pulse frequencies.
Anti scattering grid (in and out)	See Table 1a or 1b.
Parameter	Horizontal beam direction 25 cm x 25 cm x 20 cm PMMA phantom centered on top of the table / X-ray tube OT (overtable)

Table 1a dosimetry with anti-scatter grid

Mode	Settings	Resulting exposure parameters	Resulting exposure parameters	Resulting exposure parameters
General		zoom 22cm SID 900mm isocenter: 750mm interventional reference point: 600mm with 15/80 antiscatter grid	zoom 22cm SID 1050mm isocenter: 750mm interventional reference point: 600mm with 15/80 antiscatter grid	zoom 22cm SID 1200mm isocenter: 750mm interventional reference point: 600mm with 15/80 antiscatter grid
Fluoro -	kV plateau 81kV max pulse width 8ms Cu mm min: 0.6mm Cu mm max: 0.9mm Dose per frame: 32nGy 15 fps	tube voltage: 77kV tube current: 100mA exposure time: 8.0ms copper filter: 0.9mm air kerma entrance dose: 110μGy/s	tube voltage: 77kV tube current: 150mA exposure time: 8.0ms copper filter: 0.9mm air kerma entrance dose: 160μGy/s	tube voltage: 77kV tube current: 200mA exposure time: 8.0ms copper filter: 0.9mm air kerma entrance dose: 200μGy/s
Fluoro	kV plateau 70kV max pulse width 8ms Cu mm min: 0.2mm Cu mm max: 0.6mm Dose per frame: 45nGy 15 fps	tube voltage: 68kV tube current: 140mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 340μGy/s	tube voltage: 68kV tube current: 190mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 460μGy/s	tube voltage: 68kV tube current: 250mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 600μGy/s
Fluoro +	kV plateau 70kV max pulse width 8ms Cu mm min: 0.1mm Cu mm max: 0.3mm Dose per frame: 45nGy 10 fps	tube voltage: 66kV tube current: 180mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 260μGy/s	tube voltage: 66kV tube current: 240mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 340μGy/s	tube voltage: 68kV tube current: 200mA exposure time: 8.0ms copper filter: 0.2mm air kerma entrance dose: 480μGy/s
Card	kV plateau 81kV max pulse width 6.4ms Cu mm min: 0.0mm Cu mm max: 0.0mm Dose per frame: 170nGy 15 fps	tube voltage: 81kV tube current: 140mA exposure time: 6.4ms copper filter: 0.0mm air kerma entrance dose: 1700μGy/s	tube voltage: 81kV tube current: 190mA exposure time: 6.4ms copper filter: 0.0mm air kerma entrance dose: 2200μGy/s	tube voltage: 81kV tube current: 240mA exposure time: 6.4ms copper filter: 0.0mm air kerma entrance dose: 2800μGy/s
DSA	kV plateau 70kV max pulse width 160ms Cu mm min: 0.0mm Cu mm max: 0.9mm Dose per frame: 3600nGy	tube voltage: 67kV tube current: 410mA exposure time: 160ms copper filter: 0.1mm air kerma entrance dose: 3000μGy/s	tube voltage: 70kV tube current: 340mA exposure time: 160ms copper filter: 0.0mm air kerma entrance dose: 4800μGy/s	tube voltage: 70kV tube current: 440mA exposure time: 160ms copper filter: 0.0mm air kerma entrance dose: 6200μGy/s

Table 1b dosimetry without anti-scatter grid

Mode	Settings	Resulting exposure parameters	Resulting exposure parameters	Resulting exposure parameters
General		zoom 22cm SID 900mm isocenter: 750mm interventional reference point: 600mm without 15/80 antiscatter grid	zoom 22cm SID 1050mm isocenter: 750mm interventional reference point: 600mm without 15/80 antiscatter grid	zoom 22cm SID 1200mm isocenter: 750mm interventional reference point: 600mm without 15/80 antiscatter grid
Fluoro -	kV plateau 81kV max pulse width 8ms Cu mm min: 0.6mm Cu mm max: 0.9mm Dose per frame: 32nGy 15 fps	tube voltage: 77kV tube current: 100mA exposure time: 3.2ms copper filter: 0.9mm air kerma entrance dose: 45μGy/s	tube voltage: 77kV tube current: 100mA exposure time: 4.5ms copper filter: 0.9mm air kerma entrance dose: 60μGy/s	tube voltage: 77kV tube current: 100mA exposure time: 6.0ms copper filter: 0.9mm air kerma entrance dose: 80μGy/s
Fluoro	kV plateau 70kV max pulse width 8ms Cu mm min: 0.2mm Cu mm max: 0.6mm Dose per frame: 45nGy 15 fps	tube voltage: 65kV tube current: 130mA exposure time: 8ms copper filter: 0.6mm air kerma entrance dose: 100μGy/s	tube voltage: 65kV tube current: 170mA exposure time: 8ms copper filter: 0.6mm air kerma entrance dose: 130μGy/s	tube voltage: 65kV tube current: 230mA exposure time: 8ms copper filter: 0.6mm air kerma entrance dose: 180μGy/s
Fluoro +	kV plateau 70kV max pulse width 8ms Cu mm min: 0.1mm Cu mm max: 0.3mm Dose per frame: 45nGy 10 fps	tube voltage: 66kV tube current: 100mA exposure time: 4.5ms copper filter: 0.3mm air kerma entrance dose: 90μGy/s	tube voltage: 66kV tube current: 100mA exposure time: 6.5ms copper filter: 0.3mm air kerma entrance dose: 120μGy/s	tube voltage: 66kV tube current: 100mA exposure time: 8.0ms copper filter: 0.3mm air kerma entrance dose: 160μGy/s
Card	kV plateau 81kV max pulse width 6.4ms Cu mm min: 0.0mm Cu mm max: 0.0mm Dose per frame: 170nGy 15 fps	tube voltage: 81kV tube current: 90mA exposure time: 3.2ms copper filter: 0.0mm air kerma entrance dose: 600μGy/s	tube voltage: 81kV tube current: 100mA exposure time: 4.0ms copper filter: 0.0mm air kerma entrance dose: 780μGy/s	tube voltage: 81kV tube current: 100mA exposure time: 5.5ms copper filter: 0.0mm air kerma entrance dose: 1000μGy/s
DSA	kV plateau 70kV max pulse width 160ms Cu mm min: 0.0mm Cu mm max: 0.9mm Dose per frame: 3600nGy	tube voltage: 64kV tube current: 300mA exposure time: 160ms copper filter: 0.3mm air kerma entrance dose: 750μGy/s	tube voltage: 64kV tube current: 400mA exposure time: 160ms copper filter: 0.3mm air kerma entrance dose: 990μGy/s	tube voltage: 65kV tube current: 370mA exposure time: 160ms copper filter: 0.2mm air kerma entrance dose: 1500μGy/s

Table 2 Dose factor for pulse frequency



Pulse frequency [1/s]	Scaling factor for air kerma entrance dose
0.5	0.03
1.0	0.07
2.0	0.13
2.5	0.17
5.0	0.33
7.5	0.50
10.0	0.67
15.0	1.00
30.0	2.00

Table 3 Dose factor for zoom formats

Zoom format	Scaling factor for air kerma entrance dose
39	0.57
32	0.70
25	0.84
20	1.10
16	1.40
11	2.20

8. The list of components and accessories


The angiographic system Artis one			
No	Name of component	Description	Photograph
Basic block:			
1.	C-arm	The floor stand is mounted on a swivel base. It can be moved from the work position to the patient position and vice versa on the swivel base. The floor stand can be moved to different positions relative to the patient. The system allows you to obtain oblique projections using rotational and orbital motion (tilt) of the C-shaped arc.	
2.	X-ray tube assembly with primary collimator	The X-ray tube unit for cardioangiography with a primary collimator is mounted on a C-shaped arm. The collimator limits the radiation field according to the size of the region of interest (organ) or the maximum input field corresponding to the magnification level.	
3.	Flat panel detector	Image acquisition is carried out using a flat detector, followed by image processing and documentation. PD 26x30 (as30) with an input field of about 26 cm x 29 cm (39 cm diagonal) mounted on a C-arm.	


4.	Generator	X-ray generator, with a maximum output power of 100 kW.	
5.	Patient table	The Artis one system is equipped with a floor patient table. The height of the table and the longitudinal and transverse positions of the table deck can be continuously adjusted. The patient table can also be rotated. Guide rails for control modules and other accessories are attached to the table.	

6.	Mattress for patient table	The mattress is used for patient positioning, waterproof.	
7.	Cover for a mattress	Shelter for a mattress for hygienic purposes.	
8.	Fix mounted display ceiling suspension system DCS Artis one in variants: - with 30" display - with Large display 55"	Ceiling suspension for monitors in the X-ray operating room. Using the handle, you can move, rotate and adjust the height of the DCS. ARTIS one 55" Large Display* Up to three internal and four external image sources. Display layouts are configurable. Important images can be scaled to the desired size, less important information can be moved out of the focus.	 



			
9.	Footswitch for switching on radiation, in variants: - wireless with the charger; - wired	Footswitch for switching on radiation with configurable keys. A wireless 4-pedal footswitch for release of fluoroscopy, exposure, and table brake, as well as a configurable additional function.	 
10.	Control console	Modular control console, control modules are equipped with a joystick and several keys.	
11.	Imaging system unit with basic software	The imaging system consists of a computer with electronics for image creation and processing, input devices, CD / DVD drives	




			 
12.	Display in the control room	Monitors are used for medical purposes for real-time diagnosis (21-inch diagonal).	
13.	Keyboard	Computer keyboard with a block of keys and symbols for the imaging system.	



14.	Computer mouse	wheel	The mouse has two buttons and a wheel which has also a button function.	
15.	Second op. handswitch	contr.	Manual radiation switch for triggering the exposure series and single exposure.	
16.	Dose measuring chamber KermaX plus 120-160 EAS		Radiation dose measuring chamber, manufactured by IBA Dosimetry GmbH.	



17.	The operating instructions in Russian for KermaX plus 120-160 EAS	User manual for a camera for measuring radiation dose.	 KermaX® plus 120-160 EAS Руководство пользователя оборудования
18.	The operating instructions for ARTIS one in Russian	Artis one user manual in Russian.	 ARTIS one Руководство пользователя VA1 SIEMENS Healthineers



Accessories:		
1.	DR mode	X-ray mode for obtaining both series (scenes) and individual images.
2.	Card acq. mode w/high speed	Angiocardiology mode for obtaining both series (scenes) and individual images.
3.	DSA acquisition mode	In the digital subtraction angiography mode, radiographs of the examined area are performed without a contrast agent (image-mask) and with a contrast agent (images with filling).

4.	PERISTEPPING (native)	Radiography of a linear examination area (for example, a leg or an arm) can be performed in stages (without subtraction) using angiography using a stepped method.	
5.	PERISTEPPING PERIVISION	/ Step method implemented by the company	
6.	HeartSweep	Siemens through the use of digital online subtraction. When a contrast agent is injected, angiography is performed using a step method in a subtracted form.	

7.	3D Imaging & Processing	Imaging mode of the coronary arteries in different projections when performing	
8.	LA Segmentation	Left atrial segmentation mode for planning and performing electrophysiological procedures.	
9.	Mem. expansion 3 (50k - 1k Matrix)	Increase the free memory space of the master database to prevent insufficient system resources.	
10.	ECG-triggered FL	ECG synchronized fluoroscopy mode to control frame rate during fluoroscopy.	

11.	Tube Shielding for EP	Metal filters are used as shielding against radiation and to reduce the distribution of signal components in the image caused by systems electromagnetic localization, as well as ablaters. Such distortions occur only with EP interventions (electrophysiological). Made of nickel alloy. Reducing the number of bands is done with filters specially designed for this purpose. Optional color monitor for ECG data will be mounted on DCS in addition to existing one.	  
12.	CLEARstent, CLEARstent Live + CLEARstent	X-ray mode and isolation of fine structures, for example, deployable stents.	

			
13.	Lower body radiation protection	Lower body radiation shielding is designed to reduce scattered radiation in the lower body of the examining physician.	
14.	Fixed upper body rad. protection	A radiation shield is used to reduce the effects of scattered radiation on the upper body of the doctor. This is especially true of the effect of radiation on the eyes and thyroid gland of the examiner.	

			
15.	ECG radiolucent cable	ECG system connection cable for monitoring the patient's ECG or for performing X-ray or fluoroscopy synchronized with the ECG. Cable length - 30 m.	 Metal rack which holds the ECG at the bottom to fix it.
16.	VOLCANO s5i Cable Kit	Connection cables for the ultrasonic unit of the intravascular ultrasound system. Cables lengths - 30 m.	


17.	Vascular analysis	2-dimensional quantitative vessel analysis to measure stenosis in non-coronary vessels.	
18.	Scientific QCA	2-dimensional quantitative analysis of coronary vessels to measure stenosis in the coronary arteries.	

19.	QCA bifurcation	2-dimensional quantitative analysis of coronary vessels to measure stenosis of bifurcated vessel segments.	 The screenshot displays the Siemens QCA software interface. At the top, there are icons for a keyboard and a mouse. Below these, the main menu includes 'Analyze', 'Adjust', and 'Tools'. The 'Analyze' menu is open, showing three options: 'Bifurcation', 'Single', and 'Double'. The 'Bifurcation' option is selected, and a sub-menu is visible below it, showing 'Bifurcation' and 'Single'. The background of the software interface is dark blue.
-----	-----------------	--	---


20.	LV Analysis	Quantitative analysis to determine movement of the left ventricular wall hearts.	
-----	-------------	--	---


21.	IZ3D	3-dimensional quantitative analysis of coronary vessels to measure stenosis in coronary arteries in 3 dimensions from two 2-dimensional images.	
-----	------	---	---


22.	Fluoro Loop	Scene save mode of the last fluoroscopy.	
23.	Automap	The mode allows you to restore the position of the system from the selected reference image or display the reference	
24.	DICOM RIS-Modality Worklist	HIS / RIS-wide search mode by querying a worklist to access patient data destined for other workstations or modalities.	

25.	DICOM MPPS	The mode of documentation and tracking of all procedures performed.	
26.	DICOM Print	A mode for exposing images to film or printing.	
27.	StreamLink	For connecting the Artis One system to an external network.	
28.	Panoramic Display	Optional LCD display for ceiling suspension in OR.	

29.	Display connection kit	Cables for connecting and installing an additional LCD display on the ceiling suspension in the operating room. Cables lengths - 36 m.	
30.	UPS device/table/imaging system	The UPS is a battery backup system for providing power in the event of a utility power failure.	
31.	Injector+Handswitch kit for c-room	The injector is mounted on the wall using a special clip.	

32.	Injector Connection - Pedestal	Connector for connecting and installing the injector in the control room.	
33.	Intercom – Comfort Intercom ceiling LED	An intercom is used for audio communication between patient and operator when working with Artis one in the control room. Cables lengths - 20 m.	
34.	syngo Security Package (SW lic.)	To protect data at the patient and exam level.	



35.	3D resolution test	To check the quality of 3D images.	
36.	Infusion bottle holder	For installing an infusion set.	
37.	Cable clips ECG	For fastening ECG cables.	
38.	Body Module	Has guide rails for accessories and can be moved.	
39.	Head-end holder	For attaching standard accessories.	

40.	Instrument tray	The instrument table can be mounted on either side of the table on all Artis one systems. For ease of use, it can be turned and adjusted in height.	
41.	Head wedge	Used to support the head during head examinations.	
42.	Head position (prone - adult)	During invasive procedures requiring the patient to be positioned on a table in a prone position, the head can be placed on a head pillow filled with soft gel.	

43.	Head holder, deep	Used to securely support a wide head during head examinations.	
44.	Head holder, flat	Used to securely support a narrow head during head examinations.	

45.	Handle	Provides a secure fit for the patient	
46.	Handgrips support	with This accessory allows the patient to hold onto the handrails.	
47.	Arm rest	Hands rest comfortably on the supports.	
48.	Unilateral arm rest	For positioning the patient's arm without using metal parts.	

49.	Handgrip	For comfortable hand positioning. For the arm, height and length adjustment is possible. For a comfortable position of the arms along the body. Two straps with Velcro fasteners to secure the patient. They can be attached to the table top.	
50.	Arm holder	For comfortable hand positioning. For the arm, height and length adjustment is possible. For a comfortable position of the arms along the body. Two straps with Velcro fasteners to secure the patient. They can be attached to the table top.	
51.	Body strap set	For comfortable hand positioning. For the arm, height and length adjustment is possible. For a comfortable position of the arms along the body. Two straps with Velcro fasteners to secure the patient. They can be attached to the table top.	
52.	Disposable cover for model 6272	Suitable for Round Cut Radiation Shielding Screen.	
53.	Sterile cover Gr1 (99) D=50cm	Round shape, diameter 50 cm.	

54.	Sterile cover Gr2 (99) 75x80cm	Rectangular shape, 75 cm x 80 cm.	
55.	Sterile cover Gr3 (100) 90x90cm	Rectangular shape, 90 cm x 90 cm.	
56.	Sterile cover Gr4 (99) D=115cm	Round shape, diameter 115 cm.	
57.	IV-LINQ Angio- IVUS/OCT Co-Reg	Provides a joint display of the angiographic system and IVUS / OCT. The program mode comes as a download link. The user can use the existing PC for installation. In addition, a software connection tool called CAAS DICOM Connect is provided to receive DICOM files sent by Siemens Healthineers X-ray system or IVUS / OCT system.	
58.	Additional License CAAS IV-LINQ	Allows installation of IV-LINQ on multiple devices in an organization. License activation is performed automatically if the computer is connected to the Internet. If not, then the ".tmr" file that is generated when the product key is entered into the software must be sent for activation.	

59.	Video connector box	The video connector block is used to connect one video source of another type and convert those signals to DVI-D format for display on Artis displays.	
60.	Raysafe	The intend use of Raysafe i3 is designed for personal use on hospital staff. It measures scattered radiation in examination rooms using a dosimeter, and shows a graphical visualization of dose rate in real time on a display.	
61.	Spacer rail	Multi-purpose accessory rail for attachment of tableside rail.	


9. Capabilities of the Medical Device and ways of integration with other medical devices

This product may be used only in combination with software authorized by Siemens.

To ensure required operational safety, the system may be used with devices/components, application of which is explicitly authorized by Siemens.

10. Labeling information


The following system labels are placed on the bottom of the generator/system cabinet. All country specific labels necessary for e.g. registration or import in respective countries are attached to the cabinet, but are not listed here.

(1) System cabinet bottom left:

ARTIS one System identification label (example)


According to MDD:

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaixin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

According to MDR:


EC REP Siemens Healthcare GmbH, Henkestr. 127
91052 Erlangen, Germany

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaixin C. Ave., 2nd, Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Made in China

Label nearby to system identification label

Siemens
Healthcare GmbH
Henkestr. 127
91052 Erlangen
Germany

RoHS importer label

ARTIS one

COMBI

(01) 04056869009964 GTIN
(240) 10848600 MODEL
(20) 00 REV

Siemens Shenzhen Magnetic Resonance Ltd., Siemens MRI Center,
Gaixin C. Ave., 2nd Hi-Tech Industrial Park, 518057 Shenzhen,
PEOPLE'S REPUBLIC OF CHINA

Label for combination of medical devices (example)

- CONFIG: Configurable device
- MD: Medical device
- SN: Serial number
- UDI: Unique device identification
- COMBI: Combination

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
Generator connection label

Second set of labels


**SECOND SET OF LABELS
TUBE HOUSING ASSEMBLY**

Second set of labels for X-ray tube assembly on the inner side of C-arm

**SECOND SET OF LABELS
COLLIMATOR**

Second set of labels for collimator on the collimator cover

**SECOND SET OF LABELS
DETECTOR**

Second set of labels for detector on the inner side of C-arm


**SECOND SET OF LABELS
DAP CHAMBER**

Second set of labels for DAP chamber on the inner side of C-arm

MD components labels

MEGALIX Cat Plus 125/40/90-GW

pixium 2630S detector

Polydoras A100 Plus generator

Patient table


DCS for 30" display

DCS for Large display

Wire footswitch

Wireless footswitch

Wireless footswitch charger


Control units (stand, table, CCM, EHC)

MONITOR COLOR 30" 4MP DSC3012-DC


MONITOR COLOR 21.3" MX215

Large Display DSC55-BT-P


Keyboard/mouse


Handswitch

KermaX plus 120-160 EAS

MD accessories labels

Tube Shielding for EP


Lower body radiation protection


Fixed upper body rad. Protection

ECG radiolucent cable


VOLCANO s5i Cable Kit

StreamLink

Panoramic Display


Display connection kit


UPS device/table/imaging system


Injector+Handswitch kit for c-room

Injector Wall Connection – Pedestal

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















Warning label for room height

Warning label for ionizing radiation – China

This product complies with DHHS regulations
21 CFR Subchapter J, applicable at date of
manufacture.

Manufactured: Month Year
Siemens Shenzhen Magnetic Resonance Ltd.,
Siemens MRI Center
Gaoxin C. Ave., 2nd, Hi-Tech Industrial Park
518057 Shenzhen, China

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Pediatric use remove grid label – China


Plate with EAC Label*

*This symbol is applicable to the system until 11.12.2022, after this time the symbol is not applied.

Pictograms and labels

The following pictograms and labels may be present on and apply to your system.

Pictogram/ label	Short description
	Degree of protection type B (→ Page 62 Protection against electric shock)
IPX2, IPX4, IPX8	IP codes
	Equipotentiality
Duty Cycle max. 15% 90s / 10 min	Restricted duty cycle
	Power ON for part of equipment
	Power ON/OFF push button
	Maximum patient weight (→ Page 46 Maximum patient weight)

Maximum patient table load

CPR position

Non-ionizing electromagnetic radiation

Caution

Operation instructions

Consult the accompanying documents

Accessory load (example)


Maximum accessory load (example)

Collision area

Emergency STOP

Footswitch

Equipment disposal

Warning signs

Danger of crushing

This warning sign indicates possible danger of crushing for the patient and/or examiner.

Siemens Healthineers angiography systems generally possess in their standard configuration a safety system for preventing collisions (collision protection).

General danger

However, under certain circumstances, additionally installed components, like lead protective screens, lamps, auxiliary equipment etc. can cause a collision with the angiography system.

To protect these components and also to protect the patient, such components are provided with a warning sign shown on the left.

We can provide you with these warning signs. Customer Service will also attach them to additionally installed components on request.

Observe the documents

These pictograms indicate that special instructions can be found in the accompanying documents.

Please read the Operator Manual and other accompanying documents


Safety information - maximum patient weight

The permissible patient weight with which you can load your tabletop is stated on the label on the tabletop and on the table column and in the technical data.

The value on this label indicates the maximum patient load, without accessory load and without CPR load.

(CPR = Cardio-Pulmonary Resuscitation)

The value on this label indicates patient load including accessory load but without CPR load.

It is important that the load is distributed uniformly over the tabletop. Otherwise there is a risk of material deformation and system malfunction. The patient weight includes any parts permanently or loosely connected with the patient, such as equipment, prostheses, implants, plaster casts.

Example of incorrect use with uneven weight distribution:

A patient with a maximum weight sitting at the end of a fully extended tabletop.

 WARNING
Patient is too heavy for the table. Personal injury due to broken and falling parts of the table, patient falling down from the table or malfunctions of the table ◆ Check both of the labels on the table indicating the allowable weight, i.e. • Label on the table base • Label on the tabletop ◆ Do not position a patient with more weight than indicated in either of both labels.

Accessory load

Additional loads may be attached to any accessory rail or to all accessory rails in total. It includes all accessories, e.g. injector, radiation protection devices!

The weight of accessory parts is stated by labels.

(Examples)

The component, e.g. tabletop or accessory rails, must not be loaded with more than the weight indicated.

In the interest of operating safety for the patient and unit, the permissible patient weight and accessory weight must not be exceeded.

 CAUTION
Accessories are too heavy for the accessory rails and the table. Personal injury by falling parts. The weight limits for the accessories (all accessories, e.g. injector, radiation protection devices) attached to the accessory rails of the table are as follows: • Maximum 100 kg in total • Maximum 60 kg at the side rails each • Maximum 40 kg at the foot rail

Protection against electric shock

- Protection class

Class I: The system belongs to protection class I with a type B applied part according to IEC 60601-1.

- Applied part

Type B: accessory for the purpose of positioning the patient

- AP/APG

The system is not intended for use with flammable anaesthetics (Category P) or flammable anaesthetics with oxidants (Category APG) and in oxygen rich environment.

- Equipotential equalization

Systems for which equipotential equalization is recommended must only be operated in medical facilities where supplemental equipotential equalization has been installed and tested according to the specifications in DIN 57107/VDE 0107/6.81 section 5 (Germany) or the relevant local and federal regulations. Equipotential equalization connector is designed according to IEC 60601-1, Clause 16.

- Power supply

For all products that are operated within an X-ray system, the power supply must be made through a conductor or other multipole circuit breaker installed on-site.

The room installation must comply with IEC 60364-7-710 or the corresponding national regulations.

- Power outlet

The country-specific multipurpose power outlet on the patient table is not switched off and on with the system.

Only devices which accord to IEC 60601-1 must be connected.

Do not connect a multiple socket outlet or extension cord.

Do not connect life supporting devices to the power socket of the patient table. The power socket does not have the necessary medical grade.

For the ratings of the power outlet (supply voltage / over current protection), refer to the label next to the power outlet.

- Injector connector

When the injector is removed from the table, the adapter must also be removed from the connector on the table for safety reasons.

- Additional devices

If an additional device is used in examinations, e.g. ECG unit:

Before connecting the patient to an additional device, make an additional, conductive connection between a point of the equipotential equalization, e.g. the tabletop, and the device.

For approved injectors with the Siemens Healthineers central connection cable, this additional measure is not required for connection to the patient table. The central connection cable is already equipped with an additional conductor for the above purpose.

- Covers

If socket covers are damaged, they must be replaced.

In the event of defects, for example, if a covering cap has broken off: Call Customer Service.

- Opening the units

Only authorized service personnel are permitted to open the units.

Artwork of the Medical device labeling in Russian

Additionally, a Russian label is applied with the following content (label is placed on the bottom of the system cabinet).

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11.The list of International Standards and Directives applied to the Medical Device

EN 1041:2008 +A1:2013	Information supplied by the manufacturer of medical devices
ISO 10993-1:2018	Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process - Technical Corrigendum 1 (ISO 10993-1:2009/Cor 1:2010)
ISO 12052: 2017	Health informatics - Digital imaging and communication in medicine (DICOM) including workflow and data management (ISO 12052:2017)
ISO 14971:2019	Medical devices - Application of risk management to medical devices / Corrected and reprinted in 2019-12
ISO 15223-1:2021	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1 : General requirements (ISO 15223-1:2021)
IEC 60601-1-2:2020	Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Collateral Standard: Electromagnetic disturbances - Requirements and tests (IEC 60601-1-2:2014/A1:2020)
IEC 60601-1-6:2020	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability / Combines IEC 60601-1-6 (2010-01), AMD 1 (2013-10) and AMD 2 (2020-07)
IEC 62304:2015	Medical device software – Software lifecycle process (IEC 62304 EDITION 1.1 (2015-06-01) combines IEC 62304 (2006-05) and AMD 1 (2015-06)+AC:2008+A1:2015)

IEC 62366-1:2020	Medical devices – Part 1: Application of usability engineering to medical devices (IEC 62366-1(2015-02), Corrigendum 1 (2016-07) and AMD 1 (2020-06))
IEC80001-1 :2010	Application of risk management for IT-networks incorporating medical devices - Part 1: Roles, responsibilities and activities (I EC 80001-1:2010)
IEC 60601-1:2021	Medical electrical equipment – Part 1: General requirements for basic safety and essential performance / (IEC 60601-1:2005/A2:2020). +A1:2013 +AC + A12:2014 +A2:2021)
IEC 60601-1-3:2013	Medical electrical equipment - Part 1-3: General requirements for basic safety and essential performance - Collateral standard: Radiation protection in diagnostic X-ray equipment / Combines IEC 60601-1-3 (2008-01) and AMD 1 (2013-04) and AMD 2 (2021-01)
IEC 60601-2-28:2017	Medical electrical equipment - Part 2-28: Particular requirements for basic safety and essential performance of X-ray tube assemblies for medical diagnosis
IEC 60601-2-43:2019	Medical electrical equipment - Part 2-43: Particular requirements for basic safety and essential performance of X-ray equipment for interventional procedures / Combines IEC 60601-2-43 (2010-03) and AMD 1 (2017-05) and AMD 2 (2019-10)
IEC 60601-2-54: 2018	Medical electrical equipment - Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy / Combines IEC 60601-2-54 (2009-06), Corrigendum 1 (2010-03), Corrigendum 2 (2011-06) and AMD 1 (2015-04 and AMD 2 (2018-06)

12. Sterility.

The Medical device is non-sterile.

13. Disposal.

On disposing of the complete system or parts thereof, currently valid environmental legislation must be observed.

Examples of environmentally relevant components are:

- Accumulators and batteries;
- Transformers;
- Capacitors.

For details please contact your local customer service representative.

System components hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

Class A - Waste generated during/after the use of medical devices must be classified and disposed of (disposed of) as epidemiologically safe waste, close to MSW in composition, in accordance with current legislation.

All medical waste shall be collected, used, decontaminated, placed, stored, transported, listed and disposed of in accordance with the laws on public health and epidemiological welfare.

14. Packaging

Packaging consists of wood, nails, cardboard, paper and plastic. A shipping box contains a packing list indicating names, quantities, packing dates and packer's reference number.

Each package has a label containing the following information:

- manufacturer's trademark;
- product name;
- date of manufacture.

Each group box has a label according to manufacturer's working drawings, indicating the following information:


- medical device manufacturer's name;
- net and gross weight;
- box dimensions.

ackaging signs:

Handling signs

Indication of climatic conditions

15. Transportation and storage requirements and conditions.

Ambient conditions for transport and storage

The system will be stored and installed by Siemens service. The system will be fixed installed and bolted in the respective room. Therefore local transport of the user is not possible.

16. Cleaning and disinfection

WARNING

Inappropriate cleaning of the system, e.g. letting fluids enter the interior of system, or using harsh cleaning agents.

Electrical disturbances, corrosion or surface damages and correspondingly, malfunction or failure of the system

- ◆ Use only substances for cleaning and disinfection, which are recommended.
- ◆ Do not let cleaning liquids seep into the openings of the system, e.g. air openings, gaps between covers.
- ◆ Observe the following cleaning and disinfection instructions.

16.1. Cleaning agents for all LCD displays only

The following active ingredient classes can be used for LCD displays:

Active ingredient class	Tested cleaning agents and disinfectants	Additional examples
Alcohol	Ethyl alcohol 96%	Hospiset cloth Microcide liquid
Aldehyde	Melsitt Cidex	Aldasan 2000 Kohsolin Gigasept FF
Chlorine derivatives	Terralin	Quartamon Med
Disinfectants	Taski DS5001 (Diverseylever Labs) Morning Mist Surfanios Fraicheur Citron (Anios Labs)	
Guanidine derivatives	Lysoformin	
Quaternary compounds	Incidur spray, full strength	
Commercially available dishwashing liquid	Tempo	Fairy Ultra, Palmolive
Benzine	Petroleum ether	
Pyridine derivatives	Active Spray, full strength	
Water	Tap water Distilled water	

16.2. Cleaning agents for all other system parts

Some substances contained in disinfectants are known to be **hazardous to health**. Their concentration in the air must not exceed the legally defined limit.

Always follow the manufacturer's instructions when using the disinfectants.

Below suggestion is not applicable to mattresses. Solely follow the mattress manufacturer's instructions to choose suitable disinfectants and cleaning products.

The following active ingredient classes can be used:

Active ingredient class	Cleaning agents and disinfectants
Aldehyde	Incidin Perfekt
Chlorine derivatives	Clorina
Commercially available dishwashing liquid	Tempo, Fairy Ultra, Palmolive
Guanidine derivatives	Biguacid S
Quaternary compounds	Oasis Pro
Peroxides	Sanosil
Pyridine derivatives	Bacoban WB

The following products should not be used:

- All alcohol-based products
- All chlorine-releasing products
- Benzine / phenole-based products
- Products: Banicide Advanced, Wavicide, Terralin, Cidex OPA

16.3. Special instructions for all LCD displays only

- LCD displays are very sensitive to mechanical damage.

Avoid scratches, bumps, etc.

- LCD displays are very sensitive to liquids.

Longer contact with the liquids can cause discoloration or can leave calcium residue on the surface. Droplets that run down between the panel and the frame can cause damage or total failure of the panel!

1 If possible, immediately clean any droplets of liquid.

2 If a cleaning agent is sprayed directly on the display surface, make sure that any droplets that run down are wiped off with a microfiber cloth before they reach the edge of the panel.

3 If the front of the panel is soiled, clean it with a microfiber cloth and, if necessary, with a glass cleaning agent. Clean housing parts with a cleaning agent for plastics only.

16.4. Instructions for cleaning and disinfection

Sprays can penetrate the interior of equipment. It can cause damage to electronic components and the formation of flammable mixtures of air/solution vapor.

Use gloves to avoid infection.

1 Before cleaning the system, block radiation and block unit movements.

2 Use a lukewarm detergent solution and a soft cloth to remove slight contamination.

3 Clean all contaminated parts and all parts which may or have come into contact with the patient.

4 Remove contrast medium spots best of all with warm water.

5 Remove major contamination first with a cloth soaked in alcohol, then wipe off with clear water.

- 6 After using disinfectants, always wipe off with clear water.
- 7 Remove any contamination of Velcro tapes using a plastic comb.
- 8 Keep the ventilation slots of all components unobstructed.
- 9 Regularly clean the dust off all rails and joints etc.

17. Maintenance.

To ensure the safety of the patient, operating personnel and others, tests must be performed regularly to maintain the safety and proper functioning of the product.

The maintenance intervals are specified in the Maintenance Plan in the System Owner Manual.

If more frequent inspections and/or maintenance are required by national regulations, compliance with these regulations is essential.

All parts of this system in which wear could lead to a hazard must be inspected by trained technical personnel and replaced if necessary.

Maintenance work can be performed by Siemens Customer Service for a charge.

Periodic maintenance and routine checks are necessary to ensure safe and satisfactory operation of your system.

Proof of these activities is required by the authorities in some countries.

All legally required tests must be performed at the specified time intervals.

In Germany: constancy test, stipulated under §116 of the Strahlenschutzverordnung (StrlSchV).

This work may be performed only by qualified and authorized service engineers.

Qualified in this context means that the engineers have been trained accordingly or have acquired the necessary experience in practice. Authorized means that the engineers have been authorized by the operator of the system to perform maintenance work.

Prior to the first system start-up, we recommend that you name a staff member who will be responsible for ensuring that routine checks and preventive inspection and maintenance work are performed. This employee has to archive all certificates in the "System Owner's Manual" folder.

In addition to our repair service, Siemens Healthineers also offers you a complete range of services for preventive inspection and maintenance of your system. These services can be called on as required or agreed in a flexibly drafted service contract.

If you have not received a quotation from our UPTIME Services service organization, please contact the Siemens Healthineers representative responsible for your facility.

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of 10 (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and operating safety of these products.

Please contact Siemens Healthineers about these measures early enough.

This Maintenance Plan provides information about routine checks and periodic maintenance of your ARTIS one system.

Routine checks:

Normally the system operator entrusts the clinical operating personnel with the task of performing routine checks. These include daily, weekly, and monthly inspections as well as legally required checks as described in the section containing Functional and Safety Checks of the operator manual.

Safety Checks

The following checks are intended to ensure the safety of the system. Where appropriate, preventive measures must be adopted or repairs performed. The points to be checked are essentially prescribed by laws and standards.

Object or Function	Reason	What is checked	Interval (Months)
General checks			
Overall system	Safety of patients and personnel	Visual check of the system (cover panels) for damage and sharp edges	12
Cables and how cables are laid	Protection of patient and personnel against electric shocks.	Visual check of cables and how cables are laid for safety deficiencies.	12
Accessories	Safety of patients	Check for safety deficiencies	12
System radiation protection devices	Protection of patient and personnel against radiation injuries	Check the unit's radiation safety devices and any additional safety devices that may be present (lower body radiation shield, upper body radiation shield, ceiling-mounted radiation shield, etc.), to make sure they are correctly installed and free of damage.	12
Warning labels	Protection against incorrect operation and thus protection of patients, personnel and unit.	Check for presence of the required warning labels visible to the user for operation of the system.	12
User interfaces	Protection against incorrect operation and thus protection of patients, personnel and unit.	Check for legibility and operating symbols	12
Electrical Checks			
Electrical safety	Protection of patient and personnel against electric shocks.	• Ground wire resistance measurement • Unit leakage current measurement • Patient leakage current measurement of the complete system per national regulations	12
Mechanical Checks			
Floor and ceiling mounting	Protection of patients, personnel and unit	Mounting (visual check)	12
Drives	Ensure power transmission and safety functions	Check for wear, tension, ease of movement and function	12
System movements (manual)	Safety of the system latching mechanisms	Check of brake holding and end stops	12
Movable components	Protection of Patients, Personnel and Unit	Cable lead-in, movement behavior and, if applicable, brakes	12
Accessories	Safety of patients, personnel	Check for functionality and for good condition	12
User interfaces	Protection against incorrect operation and thus protection of patients, personnel and unit.	Check for completeness	12
Safety-relevant parts subject to wear	Protection of patients, personnel and unit	Check for functionality and for good condition	12

Function Check

Emergency stop devices	Preventing the first malfunction caused by application errors and collapsing patients	Switch-off of system functions after activating the emergency stop device	12
Control devices and warning indicators	Inform the operator about relevant system conditions and overload situations	The functions of the following displays: • Radiation • Fluoroscopy • X-ray tube load • Block Movements, Block Radiation • Unit movements at the collision limit	12
System movements (motorized)	Protection of patients, personnel and unit	Safety shutdown of the movements	12
Collision protection	Prevent damage to unit components	Automatic shutdown of system movements in the collision zone (e.g. ceiling, wall, floor)	12
Object or Function	Reason	What is checked	Interval (Months)
Function test	Prevent damage to unit components	Final function test of all components	12

Preventive Maintenance

The purpose of preventive maintenance is to reduce unforeseen failures to a minimum.

Thus the prerequisites for the system to meet the assured characteristics in the long term are created.

The effects of different operating conditions (full or partial load operation, temperature, size of dust particles, humidity, gases, vapors) are checked and the condition of parts subject to wear is determined by recording and analyzing characteristic values.

Preventive measures must be adopted or repairs must be undertaken as appropriate.

The specified maintenance intervals correspond to the minimum requirements.

Compliance with stricter national legislation may be necessary.

Object or Function	Reason	What is checked	Interval (Months)
Overall system	Preventive measures to avoid • Safety risks • Overheating • Wear and tear • Image artifacts	Check of operating data Check of cables, cable routing and cable connections for damage Cleaning contrast medium, blood, disinfectants from areas not accessible to the customer Inspection and cleaning of coolant and air circulation passageways Inspection and cleaning of optical transmission paths Removal of foreign objects, e.g. positioning aids and injection needles Check of the tolerance limits Check of the movement forces Check of the drive characteristics, acceleration and deceleration movements Measures to ensure that all components move smoothly Check and analysis of points of friction Repairing minor damage	24

Quality and Function Tests

Quality and function tests are used to check whether the system meets the assured characteristics. Image quality tests determine differences from the original status (e. g. resolution, contrast range, minimum contrast, image signal and, if applicable, check of Digital Subtraction Angiography).

If there are differences, preventive measures must be taken or repairs performed, whenever appropriate.

Object or Function	Reason	What is checked	Interval
Overall system	Optimum operation on the basis of the specifications listed in the data sheet	Interaction of all system components according to the assured characteristics	12 months
Vacuum components: • X-ray tube assembly	Assuring the system specifications (vacuum components are subject to aging)	Image Quality	12 months
Beam geometry, centering, beam collimation	Observance of the specifications and legal regulations to minimize radiation exposure for the patient and for personnel. (DHHS and country-specific regulations must be observed.)	Centering the central beam onto the center of the image receptor Coincidence of radiation field size and image receptor size or light field size and radiation field size	12 months
Radiation dose	Compliance with specifications and legal requirements to minimize radiation exposure of patient and personnel. (DHHS and country-specific regulations must be observed.)	Check of the dose rate/cutoff dose in • Fluoroscopy (DL) • Acquisition	12 months
Detail recognition	Compliance with specifications and legal requirements. (DHHS and country-specific regulations must be observed.) Ensuring image quality	Resolution during • Fluoroscopy (DL) • Acquisition	12 months
Image contrast	Compliance with specifications and legal requirements Ensuring image quality	Minimum contrast and dynamic range during • Fluoroscopy (DL) • Acquisition	12 months
DSA Device	Compliance with specifications and legal requirements Ensuring image quality	DSA scene: DSA scene; contrast sensitivity, dynamic range, logarithming, subtraction, artifacts	12 months
Image display	Compliance with specifications and legal requirements Ensuring image quality	Brightness and contrast of the configured monitors	12 months
Image Disturbances	Compliance with specifications and legal requirements Ensuring image quality	Image display for intolerable image artifacts for all configured application techniques	12 months
Object or Function	Reason	What is checked	Interval
Dose measuring device	Maintenance of specifications (DHHS and country-specific regulations must be observed.)	Accuracy of the display	12 months
Image documentation systems	Compliance with specifications and legal requirements. (DHHS and country-specific regulations must be observed.) Ensuring image quality	Gray scale reproduction, geometric imaging characteristics, resolution, optical density, artifacts	12 months (or in accordance with specifications from other manufacturers)

17.1. Service contacts

Regular inspections should be a part of the annual maintenance performed by Siemens Service under a service contract. We therefore recommend that you take out a service contract.

It is also possible to take out a service contract based on system usage.

Please contact your authorized local representative in this matter.

The name and address of the local representative in Russia:

Siemens Healthcare LLC
Dubinskaja str. 96
115093 Moscow, Russia
Tel.: 8 800 550 50 22
www.siemens-healthineers.com/ru/

Customer support center
8 800 550 50 22
(24 hour free line)
service.ru@siemens-healthineers.com

18. Additional service information

18.1. Pre-Installation Instructions

The medical device is not intended for installation and commissioning directly by the end customer.

Installation and commissioning of the medical device is carried out by representatives of Siemens Healthineers.

The information presented in this chapter is a partial extract from the service documentation intended for use by service specialists of Siemens Healthineers.

18.1.1. Prerequisites

Floor stand and table are each installed on a mounting plate. This chapter describes 2 different types of installation of the mounting plates:

- Installation directly on solid concrete
- Installation on filled concrete

If there is self-leveling grout or other types of floor that do not have the required load capacity as defined by the Planning Guide, the non load-bearing material must be removed in the area of the mounting plates and replace with filled concrete.

Using the heavy-duty anchor bolts included in the shipment, installation is possible in filled concrete up to a max. of 70 mm thick. Longer heavy-duty anchor bolts are available from the Hilti Company on request.

The mounting plates are secured to the floor using Hilti heavy-duty anchor bolts. The minimum required quality of the concrete is C20/25.

If needed, the mounting plates can also be installed using studs inserted through the floor into the room below. M12 or 1/2" threaded rods, minimum hardness 8.8 as defined by the ISO Norm, with large plate washers and two nuts for each on the underside of the floor are recommended. The threaded studs, plate washers and nuts must be obtained locally. Torque value of the M12 nuts: 88 Nm.

18.1.2. Preparatory work steps

Making the marks


Fig. 1: Making the marks

- 1 Floor stand mounting plate
- 2 Pilot hole
- 3 Table mounting plate
- 4 Isocenter in the head-end position
- 5 Mark A for use of the installation gauge
- 6 Control dimension

- Make the following marks on the floor in accordance with the Project Plan and Fig. 1:
 - Auxiliary line L1
 - If an installation gauge is not used: Isocenter for the floor stand in the head-end position Fig. 1 (4)
 - If an installation gauge is used:

Mark A Fig. 1 (5)

Installation gauge


Fig. 2: Installation gauge

The installation gauge is needed to align the mounting plates Fig. 2. The installation gauge is comprised of three aluminum plates that are shipped bolted together.

The installation gauge is bolted to the mounting plates for the floor stand and table using the screw with which the plate Fig. 2 (1) is attached Fig. 4.

The plate Fig. 2 (1) is intended for aligning the ceiling stand longitudinal rails with a biplane system and is not needed for the "Artis one" system.

Positioning the mounting plates

The back of the mounting plate for the floor stand can be identified by four blind holes Fig. 3 (1). Place the mounting plate down on its back side.


Fig. 3: Floor stand, back of installation plate

Assemble the two parts of the installation gauge Fig. 4 (1, 2) with screws and attach it to the mounting plates for the floor stand and table. When doing this, make sure that pilot hole C in the installation gauge Fig. 4 (4) lines up with pilot hole C in the mounting plate Fig. 1 (2).


Fig. 4: Positioning the installation gauge

1 Installation gauge, Part 1

2 Installation gauge, Part 2

3 Attach both parts of the installation gauge together using six M10 x 20 countersunk screws (wrench size: 6)

4 Pilot hole C

5 Four M10 x 20 countersunk screws (wrench size: 6) to attach the gauge to the floor stand installation plate.

6 Five M10 x 20 countersunk screws (wrench size: 6) to attach the gauge to the table installation plate.

To prevent drill dust from getting into the threads that are not used, cover the threads in both mounting plates with tape.

A cross slot Fig. 5 (A) and a long slot Fig. 5 (B) are present in the installation gauge to align the mounting plates. Align the two mounting plates with respect to the marks on the floor:

Mark A on the floor must be in the center of the cross slot.

Auxiliary Line L1 must be symmetrical to the long slot.


Fig. 5: Aligning the installation plates

Check the minimum distance at the head end between the midpoint of the floor stand installation plate and the wall:


Fig. 6: Head-end minimum distance (1)

18.1.3. Drilling the Mounting Plates

Floor stand mounting plate

Fig. 7: Drilling jig


Fig. 8: Floor stand installation plate

The drill jig included in the pre-installation kit for the floor stand Fig. 7 (1) is used together with the mounting plate as a drilling template. Here, the following should be observed:

- 1 The mounting plate is secured with seven heavy-duty expansion bolts.
- 2 It is preferable to drill holes 1-7 Fig. 8. If one hole cannot be drilled because of a concrete reinforcing bar, use the next closest hole in the mounting plate to drill a new hole.
- 3 No more than two holes may remain open between two mounting points.
- 4 If a hole is drilled incorrectly, the drill hole must be professionally filled.
- 5 Insert the drill jig into the 25 mm mounting holes Fig. 8 (1-7) in the mounting plate.
- 6 Drill holes 1-7 with a diameter of 18 mm. The depth of the drill hole depends on the type of installation:

— Installation directly on concrete:


Fig. 9: Floor stand, drill depth with installation on concrete

- 1 Floor stand mounting plate
- 2 Solid concrete

Installation on filled concrete:


Fig. 10: Floor stand, drill depth with installation on filled concrete

1 Floor stand mounting plate

2 Filled concrete

3 Self-leveling grout

4 Solid concrete

Table mounting plate

A different drill jig is required for the table. This is included in the pre-installation kit for the table. Insert the drill jig into the 20 mm mounting holes Fig. 11 (1) in the mounting plate.


Fig. 11: Table installation plate

Drill 4 holes with a diameter of 18 mm. The depth of the drill hole depends on the type of installation:

- Installation directly on concrete:


Fig. 12: Table, drilling depth with installation on concrete

1 Table installation plate

2 Solid concrete

– Installation on filled concrete:


Fig. 13: Table, drilling depth with installation on filled concrete

1 Table installation plate

2 Filled concrete

3 Self-leveling grout

4 Solid concrete

18.1.4. Inserting the Heavy-duty Anchor Bolts

- 1 If a floor covering is present, transfer the shape of the mounting plates onto the floor and remove the floor covering in this area. Place the installation gauge on the mounting plates and align it to the marks on the floor:
- 2 The plates must lie flat on the floor. If needed, remove any unevenness in the floor.
- 3 Blow out the drill holes.
- 4 Position the mounting plates again and align them to the floor marks A and L1.


Fig. 19: Mounting plates aligned

5 Place an insulation plate under each of the mounting plates for the floor stand and the table Fig. 20 (1) Fig. 21 (1).


Fig. 20: Floor stand, installing the insulation plate


Fig. 21: Table, installing the insulation plate


6 Recheck the alignment of the mounting plates.

7 Floor stand, inserting the heavy-duty expansion bolts:

CHECK WHICH HEAVY-DUTY WEDGE ANCHOR TYPE IS USED!

HILTI HSL4 ANCHOR

- red plastic ring at shaft
- HSL4 identification at shaft

Mounting Torque: 60Nm

HILTI HSL3 ANCHOR

- black plastic ring at shaft


Mounting Torque: 80Nm

Fig. 22: Types of expansion bolts

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/Installation on concrete:

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/Installation on filled concrete:


Fig. 23: Floor stand, inserting the heavy-duty expansion bolts

1a Hilti HSL-3 M12/50 or Hilti HSL-4 M12/50 heavy-duty expansion bolt with plate washer

1b Hilti HSL-3-G M12/100 or Hilti HSL-4-G M12/100 heavy-duty expansion bolt with plate washer

2 Additional plate washer

3 Insulation washer

4 Insulation sleeve

5 Mounting plate

6 Insulation plate

- Insert the heavy-duty anchor bolts Fig. 23; torque 80 Nm for Hilti HSL-3 and Hilti HSL-3-G; 60 Nm for Hilti HSL-4 and Hilti HSL-4-G.
- In addition, with installation on filled concrete:
To prevent the cover panel for the mounting plate from sitting on the threaded rods, cut off any excess higher than 19 mm Fig. 23.

8 Table, inserting the heavy-duty expansion bolts:


Fig. 24: Table, inserting the heavy-duty expansion bolts

1a Hilti HSL-3 M12/25 or Hilti HSL-4 M12/25 heavy-duty expansion bolt with plate washer

1b Hilti HSL-3 M12/100 or Hilti HSL-4 M12/100 heavy-duty expansion bolt with plate washer

2 Insulation washer

3 Insulation sleeve

4 Mounting plate

5 Insulation plate

- Insert the heavy-duty anchor bolts Fig. 24; torque 80 Nm for Hilti HSL-3 ; 60 Nm for Hilti HSL-4 .
- In addition, with installation on filled concrete:
To keep the cover panel for the mounting plate from sitting on the threaded rods, cut off any excess higher than 12 mm Fig. 24.

18.1.5. Final work steps

Caulking the mounting plates

To ensure that no liquids can penetrate between the floor and the mounting plate, caulk around the outside edge of the mounting plate Fig. 25 (Arrow) and the floor with caulking (e.g., silicone).


Fig. 25: Caulking the installation plate with silicon.

18.1.6. Floor Stand Ceiling Outlet

Installing the floor stand ceiling outlet

The cable outlet must be installed directly over the center of the floor stand mounting plate Fig. 59 (2).

The maximum admissible tolerance depends on the room height Fig. 58:

With a room height of 2700 mm to 3100 mm:

The tolerance from the center of the cable outlet to the mounting plate: 0 +100/-100mm:


Fig. 58: Room height

1 Ceiling outlet

2 Finished floor

RH Room height


Fig. 59: Floor stand, position of the cable

1 Floor Stand Mounting plate

2 Position of the Ceiling Cable Outlet

Installing the cable outlet in the drop ceiling

The following sketch shows how the cable outlet is installed in the drop ceiling; floor-ceiling distance 2700-3100 mm. When this is done, the M10 threaded studs must be shortened to an appropriate length.


Fig. 60: Cable outlet installation in the drop ceiling

1 Ceiling installation: pull forces on each screw < 1000 N

19. Service life.

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of 10 (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and operating safety of these products. Please contact Siemens Healthineers about these measures early enough.

20. Warranty.

Siemens guarantees that the Equipment has no defects caused by installation and, as of the date of risk transfer, is compliant with the approved specifications and has no defects in materials and production within 14 (fourteen) months as of the Equipment delivery date. Warranty terms and conditions for components made by the third parties and supplied along with the Equipment shall be detailed in the respective Order.

21. Contacts for reclamation.

The complaint may be submitted only for those issues which have not been the subject of the device acceptance procedure performed in accordance with the terms of the contract.

In case the system does not operate properly or does not perform commands as described in operating manual contact manufacturer representative service center in RF.

Contacts for questions and reclamations in Russia:

Siemens Healthcare LLC
Dubinskaja str. 96
115093 Moscow, Russia
Tel.: 8 800 550 50 22
www.siemens-healthineers.com/ru/

Customer support center
Tel.: 8 800 550 50 22
(24-hour free line)
service.ru@siemens-healthineers.com

Contacts for questions and reclamations in the country of legal manufacturer:

Siemens Healthcare Magnetic Resonance Ltd.
Siemens MRI Center, Gaoxin C. Ave.
2nd, Hi-Tech Industrial Park
518057 Shenzhen
PEOPLE'S REPUBLIC OF CHINA
www.siemens-healthineers.com

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This page is intended to be empty.

Notice

This OEM (Original Equipment Manufacturer) Manual is an integral part of the **KermaX-plus®**. To ensure proper use of this product, please read this manual carefully and keep it for the future reference. Do not carry out any adjustment or procedure other than those described in the manual. The attempt to do so may result in hazardous situation such as fire, explosion or electric shock to patient, operator or service engineer.

KermaX-plus® and its accessories must not be used for any other purpose than described in the accompanying documentation (intended use). Violation will result in loss of warranty.

IBA Dosimetry GmbH does not accept liability for injury to personnel or damage to equipment that may result from misuse of this equipment, failure to observe the hazard notices contained in this manual, or failure to observe local health and safety regulations.

IBA Dosimetry GmbH shall under no circumstances be liable for incidental or coincidental damage arising from use of the equipment described in this document.

No part of the accompanying documentation may be translated or reproduced without written permission of IBA Dosimetry GmbH, unless reproduction is carried out for the sole purpose to be used by several people in the same department.

The user must treat the accompanying documentation like any other copyrighted material. Especially, if part of the accompanying documentation is provided in electronic form, these files shall not be modified in any way. IBA Dosimetry GmbH and its suppliers retain title and all ownership rights to the accompanying documentation including those of translation into other languages, either in electronic or printed form.

Several key functions of the **KermaX-plus®** and its associated components are protected by international patents.

KermaX-plus® is a registered trademark of IBA Dosimetry GmbH. All other brands, trademarks and trade names mentioned in this documentation are the property of the applicable manufacturers and companies. IBA Dosimetry GmbH renounces all claims of possession to rights in brands, trademarks and trade names mentioned in this documentation, which are not their property.

IBA Dosimetry GmbH
Bahnhofstrasse 5
DE-90592 Schwarzenbruck
Germany
Phone: +49 9128 607- 0
Fax: +49 9128 607- 10
www.iba-dosimetry.com

Last update: 2021-01-15

Doc. ID: P-14-010-510-001 06 UM-KermaXplus_120-160

© 2021 IBA Dosimetry GmbH, DE-Schwarzenbruck

This page is intended to be empty.

Table of Contents

Notice	i
Table of Contents	iii
1. Introduction	1
1.1. Intended Use	1
1.2. Product Description	1
1.3. About this Manual	1
1.3.1. Conventions	2
1.3.2. Intended User	2
1.3.3. Symbols Used	3
2. Health and Safety Information	5
2.1. General	5
2.2. Important Information	5
2.3. User Responsibility	6
2.4. Regulations	6
2.5. Safety Precautions	7
2.5.1. Device Handling	8
2.5.2. Protection from Fluids	8
2.5.3. Combination with Other Systems	8
2.6. Electromagnetic Compatibility (EMC)	9
2.6.1. Electromagnetic emissions compliance	10
2.6.2. Electromagnetic Immunity compliance	11
2.6.3. Immunity to proximity fields from RF wireless equipment	14
2.7. Labels	15
2.8. Regulatory Requirements	17
3. Unpacking and Installation	19
3.1. Unpacking	19
3.2. Installation	19
4. Measurement	21
5. Quality Assurance	23
5.1. System Test (Stability Check)	23
5.1.1. Test Value	23
5.2. Calibration Check	23

5.2.1. Measurement Setup	24
5.2.2. Calibration Check Procedure	24
5.3. Test the drift of indicated values	25
5.3.1. Positive drift	25
5.3.2. Negative drift	25
5.4. Test tolerance and interval	26
5.5. Recalibration	26
6. Air Density Correction	29
6.1. Air Density Fluctuation in the Ionization Chamber	29
6.2. Correction Overview	29
7. Service and Technical Support	31
7.1. Warranty	31
7.1.1. Warranty, General	31
7.1.2. Warranty, Limitations	31
7.2. Maintenance of KermaX-plus® 120-160	32
7.3. Cleaning and Disinfection	32
7.4. Troubleshooting	33
7.5. Safety Inspection	33
7.6. Disposal and Recycling	33
7.7. Technical Support	34
7.8. Reporting Complaints	35
7.9. Returning device for repairs	35
8. Technical Specifications	37
8.1. Specifications of KermaX-plus® 120-160 EAS	37
8.2. Specifications of KermaX-plus® 120-160 MFD	38
8.3. Specifications of KermaX-plus® 120-160 LFD	39
8.4. Environmental Conditions and Requirements	40
8.5. Dimensions and Weight	41
8.5.1. KermaX-plus® 120-160 EAS / MFD	41
8.5.2. KermaX-plus® 120-160 LFD	43
8.6. Pin Assignment for the RJ45 8 pin Plug	45
8.7. CAN Node ID Configuration	45
8.8. Typical Energy Dependence	46
8.8.1. KermaX-plus® 120-160 EAS	46
8.8.2. KermaX-plus® 120-160 LFD / MFD	47
9. Index	49

1. Introduction

1.1. Intended Use

The **KermaX-plus® 120-160** is designed to be installed in diagnostic X-ray units (attached at the collimator) and is intended to be used for the measurement of Dose Area Product (DAP), DAP rate, and exposure time for standard and pediatric procedures and longtime fluoro applications.

The system and its accessories must not be used for any other purpose than described in the accompanying documentation (intended use).

1.2. Product Description

The **KermaX-plus® 120-160** is a round none transparent ionization chamber with integrated electronics to be used in one measuring mode to determine Dose Area Product (DAP), DAP rate, and exposure time. The integration into the system is customer specific. The chamber is attached to the output opening of the X-ray collimator.

The connection is through a RJ45 Plug directly to the imaging system of the X-ray system via CAN bus.

The communication interface, which provide DAP measurement data, is based on CANopen. The **KermaX-plus® 120-160** device as slave, server or consumer component, depending on the used communication object, provides function or service to the Master, Client or Producer which initiates requests for such services and is able to receive all measured data from the **KermaX-plus® 120-160** and is able to control the device. The **KermaX-plus® 120-160** system is answering all request with certain response telegrams. If it is not able to "understand" the command, it is sending a failure response telegram.

1.3. About this Manual

This manual contains information necessary to use the **KermaX-plus® 120-160** DAP Measurement system. It's describes safety information, functionalities and applications of the **KermaX-plus® 120-160**.

To ensure proper use of this product, please read this manual carefully and keep it for future reference. Do not carry out any adjustment or procedure other than those described in the manual. The attempt to do so may result in hazardous situation such as fire, explosion or electric shock to the patient, operator or service engineer.

The operator must be trained in the proper operation of the product.

NOTICE

IMPORTANT NOTICE

PICTURES AND SCREENSHOTS

All numbers and selections displayed in pictures and screenshots are only examples and no recommendations for settings or entries.

1.3.1. Conventions

The functions of the device, dialog captions and dialog text are indicated by **bold** font.

Examples: **Settings**, **Field size**.

Referrals to chapter and section headings in this manual are indicated by *italic* font.

Examples: *Notice*, *Technical Specifications*.

Throughout this OEM manual, hazardous situations or operations are identified by DANGER, WARNING, CAUTION and NOTICE. They are indicated by specific signs and colours, described below:

Sign	Meaning
	DANGER indicates a hazardous situation, which, if not avoided, <u>will result in death or serious injury</u> of the operator or patient.
	WARNING indicates a hazardous situation, which, if not avoided, <u>could result in death or serious injury</u> of the operator or patient.
	CAUTION, used with a safety alert symbol, indicates a hazardous situation, which, if not avoided, <u>could result in minor or moderate injury</u> of the operator or patient.
	CAUTION, without the safety alert symbol, used to address issues related to possible hardware damage.
	IMPORTANT NOTICE used to address operational issues not related to personal injury or hardware damages.

1.3.2. Intended User

This manual is intended for personnel with the following expertise:

Area	Expertise
Installation	Experts
Start-up, operation and shutdown	Experts or trained personnel
Maintenance	Experts
Troubleshooting	Experts

Trained personnel:

Any personnel who have received a training for using this medical device as described in this User's Guide by the expert or other trained personnel.

Experts:

Assigned person by IBA Dosimetry who received a specific training for this medical device as described in this User's Guide.

1.3.3. Symbols Used

The following symbols are used on the product and in this manual:

	Consult the User's Guide before use
	The device meets the essential requirements of Council Directive 93/42/EEC concerning medical devices.
	UL Recognized Component, E352291. Investigated to: • ANSI/AAMI ES60601-1 (2012) • CAN/CSA-C22.2 No. 60601-1 (2014) • IEC 60601-1 (2012)
	Symbol for Class II Equipment
	DC current
	Energy dependence, rated range of X-ray source voltage.
	Radiation filter, quality equivalent filtration
	The users are not allowed to change the original factory values of the parameters in the setup menus of the displays with this warning sign
	Recycling
	Separate collection of electrical and electronic devices in accordance with EC Directive 2012/19/EC: Do not dispose of the device with normal domestic waste. Keep separate from domestic waste and dispose in an environmentally safe way in compliance with local regulations.
	Manufacturer: this symbol shall be accompanied by the name and the address (incl. country) of the manufacturer
	Manufacturing date: this symbol is accompanied by a date to indicate the date of manufacture, as four digits for the year
	Serial Number
	Catalogue Number
	Medical Device

This page is intended to be empty.

2. Health and Safety Information

2.1. General

The purpose of this chapter is to identify the hazards associated with the equipment. This information is presented by displaying all safety and rating labels, which are attached to the equipment, and by providing instructions to avoid the associated hazards.

To ensure proper use of this product, please read this manual carefully and keep it for future reference. Do not carry out any adjustment or procedure other than those described in the manual. The attempt to do so may result in hazardous situation such as fire; explosion or electric shock to the patient, operator or service engineer. The operator must be trained in the proper operation of the product.

NOTICE

IMPORTANT NOTICE

ALL PERSONNEL MUST READ THIS CHAPTER

All personnel must read this chapter and be fully aware of its contents before commencing installation work, and before operating or servicing the **KermaX-plus® 120-160**.

If the **KermaX-plus® 120-160** is used in a way not specified in this OEM manual, the protection provided by the equipment may be reduced.

2.2. Important Information

This manual is part of the measurement system for the Dose Area Product), **KermaX-plus® 120-160** and is needed to ensure the correct functioning of this device. The **KermaX-plus® 120-160** including all the accessories must not be used for any other purpose as described in this manual.

CAUTION

CAUTION

WITHDRAWN FROM SERVICE IF A INVOLVING HAZARDOUS INCIDENT

If this measuring unit is involved or associated directly or indirectly with a hazardous incident, it needs to be withdrawn from service instantly. Report the details to your dealer or the manufacturer immediately.

NOTICE

IMPORTANT NOTICE

PATENTS

Several features of the **KermaX-plus® 120-160** and its components are subject to filed patent applications.

NOTICE

IMPORTANT NOTICE

DO NOT MODIFY THIS PRODUCT

Modifications to this product or the content of this manual will invalidate the product warranty and may compromise a safe operation. Therefore, do not under any circumstances, make any changes to the product or its components.

2.3. User Responsibility

The personnel, using the the **KermaX-plus® 120-160** to record the Dose Area Product, carry the full responsibility to judge each measurement result critically.

The **KermaX-plus® 120-160** should be stored in a clean, dry environment. Protect it from mechanical and thermal stress, dust and unnecessary moisture.

CAUTION

CAUTION

WET DEVICES

Devices on which moisture (condensation) has developed as a result of temperature changes must not be used unless they have been completely dried.

2.4. Regulations

The installer and operator are responsible for complying with all local regulations regarding installation and operation of the **KermaX-plus® 120-160**.

Please be advised that other mobile electronic devices, e.g. cellular telephones, exceeding the established emission limits in the EMC standard may disrupt the function of the device.

CE mark:

This medical device is labeled with a CE mark and complies to the 93/42/EEC directive. Its design guarantees provision to comply with all recommended standards, as with IEC 60601-1-2 for electromagnetic compatibility (EMC), on the level of the system component. The overall safeguard of EMC test, radiation protection and other safety features of the X-ray system and the respective test results have to be guaranteed and documented by the system manufacturer.

National regulations:

In all countries, the legally established regulations are to be observed.

Legally required tests:

All legally required tests must be performed at the prescribed time intervals, e.g. constancy test according to the X-ray ordinance (§16 RöV) in the Federal Republic of Germany, e.g. tests based on DHHS guidelines (Department of Health and Human Services) where applicable.

2.5. Safety Precautions

NOTICE

IMPORTANT NOTICE

LIABILITY

As manufacturer, IBA Dosimetry GmbH will not be held responsible for the safety features, reliability and performance of the system if:

- The system is used in a manner other than that specified in the operating manual.
- Installation, upgrades, resetting or repairs are performed by unauthorized personnel.
- Components affecting product safety are not replaced with original IBA Dosimetry GmbH spare parts.

- The **KermaX-plus® 120-160** should only be used by personnel which fulfils the following criteria:
 - Understanding of the functions and the limitations of the device as they relate to the measurement of radiation output.
 - Knowledgeable about safety procedures to be observed when working with radiation sources.
 - Competence of the security requirements of radiation and the relevant standards like IEC 60601-1 and others.
 - Experience with the use of measuring systems for the Dose Area Product
- Prior to the use of the measuring system, the correct mounting and connection of the chamber has to be ensured. The user has to check the general functionality and security of the chamber as well as all the connection cables (also for mechanical damages).
- The **KermaX-plus® 120-160** has to be used in a fire secured, clean, dry, and climatic room. Protect it from mechanical and thermal pressure as well as dust and unneeded humidity.

CAUTION

CAUTION

USING ONLY DRY CHAMBERS

If condensation is build up in the chamber because of temperature changes, the system has to be dried completely and should not be used.

- The **KermaX-plus® 120-160** is not designed for patient contact in its intended use.
- Do not remove any signs or labels of the **KermaX-plus® 120-160**. If labels become unreadable, the corresponding parts have to be sent back to the manufacturer, in order to be identified and to be relabelled.
- Never open the chamber device.

⚠ WARNING

WARNING

HIGH VOLTAGE

Because of the high voltage, there is a risk of life!

- The **KermaX-plus® 120-160** is not water proof. If there is a risk of spray, device must be covered and protected sufficiently.
- The chamber is not completely sealed, due to the special mounting situation. When the device is mounted and then used, make sure, that no foreign substances get into the measuring chamber. Any foreign subjects could lead to errors.

- Because of the high sensitivity of the measuring chamber, vibrations or touching the chamber could cause partial discharge, which could lead to a wrong measuring result.

This is the reason, why it has to be guaranteed in the C-shaped system, that those impulses are only evaluated during the exposure

2.5.1. Device Handling

- Before using the system, the operator must ensure that the chamber and all accessories are in proper working condition. The user must verify the general functionality, safety, and condition of the device and cables.
- Neither the **KermaX-plus® 120-160** chamber nor any peripheral device must be used in direct contact to the patient.
- Do not remove any labels from parts of the system components or its accessories.
- Do not open the chamber device in order to avoid contact with hazardous high voltage.
- The housing of the chamber device is not waterproof. When the risk of splashing fluid is present, mount it under auxiliary covers.
- The **KermaX-plus® 120-160** requires a voltage in the range of 15 - 24 V DC $\pm 20\%$, which is normally provided by a power supply with 4kV isolation according to IEC 60601-1. Connecting the chamber to another external power supply shall only be done with the approval and under the sole responsibility of the manufacturer of the diagnostic X-ray machine. The range and type of the voltage provided shall be checked and, if necessary, discussed with IBA Dosimetry GmbH.

WARNING

WARNING

ENSURE A SECURE INSTALLATION

When mounting the chamber and its accessories (holder, filter, etc.) on a X-ray machine, make sure the devices are locked securely in place to prevent falling down of the device when the collimator rotates. It may damage the device or cause personal injuries.

2.5.2. Protection from Fluids

Do not allow fluids to enter the converter either during normal operation or during cleaning and disinfection as this may damage the system or cause a system malfunction.

2.5.3. Combination with Other Systems

The use of accessory equipment not complying with the equivalent safety requirements of this equipment may lead to a reduced level of safety of the resulting system. Consideration relating to the choice shall include:

- Use of the accessory in the patient vicinity
- Evidence that the safety certification of the accessory has been performed in accordance with the appropriate IEC 60601-1 harmonized national standards

2.6. Electromagnetic Compatibility (EMC)

The electromagnetic environment of intended use is the professional healthcare facility environment, except for the RF shielded room of a medical electrical system for magnetic resonance imaging, where the intensity of electromagnetic disturbances are high.

(See chapter *Intended Use*).

Due to electromagnetic disturbances the measurement performance of dose area product meter can be degraded or lost.

NOTICE

IMPORTANT NOTICE

OPERATING IN NEAR OF HF SURGICAL EQUIPMENT

The Operation of **KermaX-plus®** equipment in a HF-SURGICAL ENVIRONMENT is generally excluded. Its use in such Environment the manufacturer of the diagnostic X-ray machine carries the full and entire responsibility to evaluate the correct operating and measurement of the **KermaX-plus®** equipment.

NOTICE

IMPORTANT NOTICE

OTHER EQUIPMENT CLOSE TO THE OPERATING DEVICE

The equipment or system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or system should be observed to verify normal operation in the configuration in which it will be used.

NOTICE

IMPORTANT NOTICE

USE OF ACCESSORIES, TRANSDUCERS AND CABLES

The use of accessories, transducers and cables other than those specified or provided by the manufacturer of the **KermaX-plus®** could result in increased electromagnetic emissions or decreased electromagnetic immunity of this equipment and result in improper operation.

NOTICE

IMPORTANT NOTICE

USE OF PORTABLE RF COMMUNICATIONS EQUIPMENT

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the **KermaX-plus®**, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

2.6.1. Electromagnetic emissions compliance

manufacturer's declaration – electromagnetic emissions		
Phenomenon	Basic EMC standard or test method	EMISSIONS class and group
radiated RF EMISSIONS <1GHz	CISPR 11 ^{b)}	Class B Group 1
Harmonic distortion	IEC 61000-3-2 ^{a)}	not applicable, no connection to public mains network.
Voltage fluctuations and flicker	IEC 61000-3-3 ^{a)}	not applicable, no connection to public mains network.
^{a)} This test is not applicable in this environment unless the ME EQUIPMENT and ME SYSTEMS used there will be connected to the PUBLIC MAINS NETWORK and the power input is otherwise within the scope of the Basic EMC standard. ^{b)} Mains conducted EMISSION not tested; direct connection to PUBLIC MAINS NETWORK not possible		

2.6.2. Electromagnetic Immunity compliance

manufacturer's declaration – electromagnetic immunity ENCLOSURE PORT		
Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
ELECTROSTATIC DISCHARGE ^{a)}	IEC 61000-4-2	not applicable, no parts accessible during normal operation
Radiated RF EM fields ^{b)}	IEC 61000-4-3	3 V/m ^{g)} 80 MHz – 2,7 GHz ^{c)} 80 % AM at 1 kHz ^{d)}
Proximity fields from RF wireless communications equipment	IEC 61000-4-3	see chapter 2.6.3
RATED power frequency magnetic fields ^{e) f)}	IEC 61000-4-8	30 A/m ^{h)} 50 Hz or 60 Hz
^{a)} This test applies to various external parts of the complete dose area measurement equipment which may be touched by the OPERATOR during a normal measurement (i.e. not to those parts of the IONIZATION CHAMBER and MEASURING ASSEMBLY that are normally exposed in the radiation beam). ^{b)} The interface between the PATIENT physiological signal simulation, if used, and the ME EQUIPMENT or ME SYSTEM shall be located within 0,1 m of the vertical plane of the uniform field area in one orientation of the ME EQUIPMENT or ME SYSTEM. ^{c)} ME EQUIPMENT and ME SYSTEMS that intentionally receive RF electromagnetic energy for the purpose of their operation shall be tested at the frequency of reception. Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. This test assesses the BASIC SAFETY and ESSENTIAL PERFORMANCE of an intentional receiver when an ambient signal is in the passband. It is understood that the receiver might not achieve normal reception during the test. ^{d)} Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS. ^{e)} Applies only to ME EQUIPMENT and ME SYSTEMS with magnetically sensitive components or circuitry. ^{f)} During the test, the ME EQUIPMENT or ME SYSTEM may be powered at any NOMINAL input voltage, but with the same frequency as the test signal. ^{g)} Before modulation is applied. ^{h)} This test level assumes a minimum distance between the ME EQUIPMENT or ME SYSTEM and sources of power frequency magnetic field of at least 15 cm. If the RISK ANALYSIS shows that the ME EQUIPMENT or ME SYSTEM will be used closer than 15 cm to sources of power frequency magnetic field, the IMMUNITY TEST LEVEL shall be adjusted as appropriate for the minimum expected distance.		

manufacturer's declaration – electromagnetic immunity
Input d.c. power PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
Electrical fast transients / bursts ^{a) g) k)}	IEC 61000-4-4	± 2 kV 100 kHz repetition frequency
Surges Line-to-line ^{a) b) g)}	IEC 61000-4-5	not applicable
Surges Line-to-ground ^{a) b) g)}	IEC 61000-4-5	not applicable
Conducted disturbances induced by RF fields ^{a) c) d) i) k)}	IEC 61000-4-6	3 V ^{h)} 0,15 MHz – 80 MHz 6 V ^{h)} in ISM bands between 0,15 MHz – 80 MHz ^{j)} 80 % AM at 1 kHz ^{e)}
Electrical transient conduction along supply lines ^{f)}	ISO 7637-2	not applicable

- ^{a)} The test is applicable to all d.c. power PORTS intended to be connected permanently to cables longer than 3 m.
- ^{b)} All ME EQUIPMENT and ME SYSTEM cables shall be attached during the test.
- ^{c)} INTERNALLY POWERED ME EQUIPMENT is exempt from this test if it cannot be used during battery charging, is of less than 0,4 m maximum dimension including the maximum length of all cables specified and has no connection to earth, telecommunications systems, any other equipment or a PATIENT.
- ^{d)} The test may be performed with the ME EQUIPMENT or ME SYSTEM powered at any one of its NOMINAL input voltages.
- ^{e)} Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.
- ^{f)} For ME EQUIPMENT and ME SYSTEMS intended to be installed in passenger cars and light commercial vehicles including ambulances fitted with 12 V electrical systems or commercial vehicles including ambulances fitted with 24 V electrical systems.
- ^{g)} Direct coupling shall be used.
- ^{h)} r.m.s., before modulation is applied.
- ⁱ⁾ If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.
- ^{j)} The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz; 3,5 MHz to 4,0 MHz; 5,3 MHz to 5,4 MHz; 7 MHz to 7,3 MHz; 10,1 MHz to 10,15 MHz; 14 MHz to 14,2 MHz; 18,07 MHz to 18,17 MHz; 21,0 MHz to 21,4 MHz; 24,89 MHz to 24,99 MHz; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.
- ^{k)} Not applicable only for KermaX-plus® 120-160 LFD.

manufacturer's declaration – electromagnetic immunity
Signal input/output parts PORT

Phenomenon	Basic EMC standard or test method	IMMUNITY TEST LEVELS
ELECTROSTATIC DISCHARGE ^{e)}	IEC 61000-4-2	not applicable, connector is not accessible during normal operation
Electrical fast transients / bursts ^{b) f) j)}	IEC 61000-4-4	± 1 kV 100 kHz repetition frequency
Surges Line-to-ground ^{a)}	IEC 61000-4-5	not applicable
Conducted disturbances induced by RF fields ^{b) d) g) j)}	IEC 61000-4-6	3 V ^{h)} 0,15 MHz – 80 MHz 6 V ^{h)} in ISM bands between 0,15 MHz – 80 MHz ⁱ⁾ 80 % AM at 1 kHz ^{c)}

^{a)} This test applies only to output lines intended to connect directly to outdoor cables.

^{b)} SIP/SOPS whose maximum cable length is less than 3 m in length are excluded.

^{c)} Testing may be performed at other modulation frequencies identified by the RISK MANAGEMENT PROCESS.

^{d)} Calibration for current injection clamps shall be performed in a 150 Ω system.

^{e)} Connectors shall be tested per 8.3.2 and Table 4 of IEC 61000-4-2:2008. For insulated connector shells, perform air discharge testing to the connector shell and the pins using the rounded tip finger of the ESD generator, with the exception that the only connector pins that are tested are those that can be contacted or touched, under conditions of INTENDED USE, by the standard test finger shown in Figure 6 of the general standard, applied in a bent or straight position.

^{f)} Capacitive coupling shall be used.

^{g)} If the frequency stepping skips over an ISM or amateur radio band, as applicable, an additional test frequency shall be used in the ISM or amateur radio band. This applies to each ISM and amateur radio band within the specified frequency range.

^{h)} r.m.s., before modulation is applied.

ⁱ⁾ The ISM (industrial, scientific and medical) bands between 0,15 MHz and 80 MHz are 6,765 MHz to 6,795 MHz; 13,553 MHz to 13,567 MHz; 26,957 MHz to 27,283 MHz and 40,66 MHz to 40,70 MHz. The amateur radio bands between 0,15 MHz and 80 MHz are 1,8 MHz to 2,0 MHz; 3,5 MHz to 4,0 MHz; 5,3 MHz to 5,4 MHz; 7 MHz to 7,3 MHz; 10,1 MHz to 10,15 MHz; 14 MHz to 14,2 MHz; 18,07 MHz to 18,17 MHz; 21,0 MHz to 21,4 MHz; 24,89 MHz to 24,99 MHz; 28,0 MHz to 29,7 MHz and 50,0 MHz to 54,0 MHz.

^{j)} Not applicable only for KermaX-plus® 120-160 LFD.

2.6.3. Immunity to proximity fields from RF wireless equipment

manufacturer's declaration – IMMUNITY to RF wireless communications equipment						
Test frequency [MHz]	Band ^{a)} [MHz]	Service ^{a)}	Modulation ^{b)}	Maximum power [W]	Distance [m]	IMMUNITY TEST LEVEL [V/m]
385	380 – 390	TETRA 400	Pulse modulation ^{b)} 217 Hz	1,8	0,3	27
450	430 – 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 – 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745						
780						
810	800 – 960	GSM 800/900, TETRA 800, iDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 217 Hz	2	0,3	28
870						
930						
1 720	1 700 – 1 990	GSM 1800, CDMA 1900, GSM 1900, DECT, LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 – 2 570	Bluetooth, WLAN 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5 240	5 100 – 5 800	WLAN 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5 500						
5 785						
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the ME EQUIPMENT or ME SYSTEM may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
^{a)} For some services, only the uplink frequencies are included.						
^{b)} The carrier shall be modulated using a 50 % duty cycle square wave signal.						
^{c)} As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

2.7. Labels

The labels are located on top and bottom of the chamber. One label indicates the identification and electrical ratings. Another label indicates the device model number and serial number. They are located on the rear side of the chamber.

An example is shown as below:

Figure 1 KermaX-plus® 120-160 variant chamber top view


Figure 2 KermaX-plus® 120-160 variant chamber bottom view

DHHS label

SIEMENS Device label

2.8. Regulatory Requirements

The **KermaX-plus® 120-160** fulfils the requirements of the Medical Device Directive 93/42/EEC. It is a medical device class IIb according to annex IX, classification rule 10.

The quality management system in IBA Dosimetry GmbH is certified according to EN ISO 13485.

The following standards apply:

- IEC 60580:2019
- UL 2601-1 and UL 60601-1 ed.1
- CSA C22.2 No. 601.1
- IEC 60601-1:2005+A1:2012
- IEC 60601-1-2:2014
- ANSI/AAMI ES 60601-1: A1:2012, C1:2009/(R)2012 and A2:2010/(R)2012
- CSA CAN/CSA-C22.2 No. 60601-1:14

The **KermaX-plus® 120-160** is manufactured by:

IBA Dosimetry GmbH
Bahnhofstrasse 5
DE-90592 Schwarzenbruck
Germany

NOTICE

IMPORTANT NOTICE

UL RECOGNIZED COMPONENT MARK

The KermaX-plus® is marked as UL Recognized Component and the Model Types are listed in the UL File Number E194025 and E352291.

NOTICE

IMPORTANT NOTICE

COMPLY WITH ALL LOCAL REGULATIONS

The installer and operator are responsible for complying with all local regulations regarding installation and operation of the KermaX-plus® measurement system with its accessories.

The overall safeguard of EMC test, radiation protection and other safety features of the X-ray system and the respective test results have to be guaranteed and documented by the system manufacturer.

RoHS and REACH Compliance:

The products specified have been designed and manufactured in compliance with the following specification.

- Compliant with the European **RoHS** Directive (**R**estriction of **H**azardous **S**ubstances), EC Directive 2011/65/EU
- Compliant with the European **REACH** (**R**egistration, **E**valuation, **A**uthorization of **C**hemicals), EC No. 1907/2006
- Compliance with the specification has been verified by internal design controls. We reserve to ourselves making use of the exception for medical devices as defined in EC Directive 2012/19/EC : (**WEEE** = **W**aste **E**lectrical and **E**lectronic **E**quipment)

3. Unpacking and Installation

3.1. Unpacking

NOTICE

IMPORTANT NOTICE

PERSONNEL QUALIFICATIONS

Only qualified personnel is permitted to unpack, install and operate the KermaX-plus® 120-160.

NOTICE

IMPORTANT NOTICE

TRANSPORTATION

During transportation the KermaX-plus® 120-160 should be adequately protected in the originally supplied or equivalent packing material.

3.2. Installation

Please refer first to the respective instruction manual from the X-ray system manufacturer.

CAUTION

CAUTION

RESPONSIBILITY OF INSTALLATION

The KermaX-plus® 120-160 is integrated Part of X-ray system, the system manufacturer carries the full and entire responsibility for correct and safety installation.

general mounting procedure:

- put the chamber on the collimator aperture and mount it's on the existing mounting fixture
- connect the RJ45 CAN bus plug of the KermaX-plus® 120-160 device with the X-ray system

Figure 4 Collimator aperture SIEMENS Artis zee with Model type EAS [top]

Figure 3 Collimator aperture SIEMENS Artis zee with Model type EAS [side]

CAUTION

CAUTIONS

SAFETY PRECAUTIONS

When installing the device, please observe the safety precautions described in Section 2.5.

CAUTION

CAUTION

DO NOT EXCEED THE ACTIVE AREA OF THE CHAMBER

The size of the radiation field at the surface of the centred chamber may not exceed the active area of the chamber (see Chapter 8, *Technical Specifications*).

NOTICE

IMPORTANT NOTICE

REFERENCE DIRECTION OF INCIDENT RADATION

The central beam direction of incident radiation must be orthogonal to the chamber surface.

4. Measurement

1. Before starting a measurement, check the System and Error State to be ready for measurement.
(see also Section 5.1)
2. Reset the last measurement values for new a session resp. study.
3. Make an exposure and read out the measured values with corresponding command.

⚠ CAUTION

CAUTION

DO NOT EXCEED THE ACTIVE AREA OF THE CHAMBER!

The size of the radiation field at the surface of the centred chamber may not exceed active area of the chamber (see Chapter 8, Technical Specifications).

⚠ CAUTION

CAUTION

MAXIMUM DOSE SURFACE PRODUCT RATE

The maximum dose surface product rate of 3000.00 $\mu\text{Gym}^2/\text{s}$ must not be exceeded.

⚠ CAUTION

CAUTION

RECALIBRATION / VERIFICATION

If a compensating filter is used behind the measuring chamber and before the patient the indicated level must be adjusted accordingly. The absorption of the table top when available should be compensated if the chamber is installed below the examination table.

NOTICE

IMPORTANT NOTICE

DAP CORRECTION FACTOR PROCEDURE

The measurement procedure for determination of DAP Correction Factor is described in Chapter 5.5 *Recalibration*.

This page is intended to be empty.

5. Quality Assurance

5.1. System Test (Stability Check)

To check the functionality and stability of the DAP-measurement system, a test signal can be activated by software command from the X-ray unit. The activation of this test in the DAP measurement system will only happen if the high voltage inside the chamber is within the correct range and no device error occurs.

5.1.1. Test Value

The test command activates a test charge in the amplifier of the connected chamber. Afterwards the chamber returns a value of 10000 ($\pm 100.00 \mu\text{Gym}^2$) $\pm 20\%$ which should be evaluated through the user (resp. X-Ray unit). The stability test should be performed before each measurement session.

5.2. Calibration Check

NOTICE

IMPORTANT NOTICE

FIRST CHECK

The DAP measuring system is calibrated by the manufacturer prior to shipment and a copy of the Calibration Certificate is included with the documentation. As a rule, the first check of the calibration normally takes place at the time of installation.

CAUTION

CAUTION

ABSORBER CORRECTION

The DAP system is calibrated without an absorber. Consider each correction that could be necessary to take into account the local conditions.

If a compensating filter is used (e.g. behind the measuring chamber or in a cosmetic cover of the X-ray tube), the indicated level must be adjusted accordingly (see Chapter 5.5 *Recalibration*).

The absorption of the table top when available should be compensated if the chamber is installed below the examination table.

5.2.1. Measurement Setup

To ensure correct measurements, the operator should check the systems in compliance with the calibration conditions in the table below.

Reference Conditions for Calibration	
Air humidity	50 % (relative humidity, no condensation)
Air pressure	1013 hPa
Field size	5 x 5 cm at chamber level (max. 9 x 9 cm)
Stabilization time	10 minutes
Temperature	20 °C
Radiation Quality	100 kV (RQR 8, IEC 61267)

5.2.2. Calibration Check Procedure

1. Perform a stability check as described in Chapter 5.1 *System Test (Stability Check)*.
2. Irradiated field size parameters should be set at the control console resp. collimator device. It is recommended to collimate an average size of the irradiated field of the active chamber surface.
3. Determine the reference field area. It is the size of the surface area on the reference point level (e.g. table top, grid assembly) which is orthogonal to the central beam direction. This is the location of the reference dosimeter.

For a more accurate field size determination, it is recommended to exposure an X-ray film or a fluoroscopic screen.

If an iris shutters / collimators are used, the area formula for regular polygons shall be applied.
4. Put the reference dosimeter on the reference point level.
5. Exposure parameters should be set at the control console and an X-ray exposure released. In order to make a valid comparison with the manufacturer's calibration (see Table above, reference conditions), the radiograph must be made using 100 kV.
6. After the exposure has been released, read and take note of the results displayed on both measuring devices.

e.g. reference dosimeter display: 2703 µGy
KermaX-plus® 120-160: 33.90 µGym²

7. Comparison of the reference dose area product (*reference dose x reference field area*) and KermaX-plus® 120-160 measured dose area product:

$$\frac{100\% \times 33.90 \mu\text{Gym}^2}{2703 \mu\text{Gy} \times 0.012 \text{m}^2} - 100\% = 4.51\%$$

KermaX-plus® 120-160 value about 4.51% higher is still within tolerance. The manufacturer's tolerance is $\pm 5\%$.

NOTICE

IMPORTANT NOTICE

AIR DENSITY CORRECTION

When using a Ion chamber as a reference system, then the reference value need an air density correction (section 6 *Air Density Correction*).

5.3. Test the drift of indicated values

5.3.1. Positive drift

During absence of radiation, and after resetting the DAP meter, the indicated values should be not more than 1 count for at least 1 hour.

After the stabilization time is elapsed and a new background compensation is completed (activated by software command), check this by reading after 15 min, 30 min, 45 min and 1 h after the DAP meter was reset, and with no resetting or compensation adjustment during the test.

NOTICE

IMPORTANT NOTICE

INFLUENCE OF ENVIROMENTAL CONDITIONS

Make sure that there is no moisture in chamber cavity, the operating conditions are met and no vibrations occur (e.g. movement of system).

5.3.2. Negative drift

Perform the system test (Section 5.1, *System Test (Stability Check)*) and repeat this test measurements separated in 5min intervals, the indicated values should not more deviate than 2% from each other.

After the stabilization time is elapsed and a new background compensation is completed (activated by software command), check this by reading after 5 min, 10 min, 15 min, 20min after the stability check of the DAP meter was initialized.

5.4. Test tolerance and interval

The radiological and electrical calibration of the DAP-measuring system is already done by IBA Dosimetry.

The check of the calibration can be done on-site by a competent person or in our calibration laboratory at IBA Dosimetry. The calibration value is saved in the chamber

The following tests should be done in regular intervals to ensure the functionality of the device.

Suggested time table and tolerance for tests:

Type of test	Frequency	Tolerance
Evaluation and routine testing	maintenance and essential modification of X-ray unit	$\pm 20 \%$
Stability check (System test)	Monthly	$\pm 20\%$
Check of calibration	• at installation • every 2 years • any case following a repair	$\pm 5\%$
Drift of indicated values	Maintenance	1 count @ positive Drift Test 2% @ negative Drift Test

NOTICE

IMPORTANT NOTICE

MAINTENANCE CHECKS

The maintenance checks specified here are only recommendations. Please review and observe national and international laws for the country of use together with the local safety regulations.

5.5. Recalibration

NOTICE

IMPORTANT NOTICE

RECALIBRATION

The change of calibration should only be carried out in exceptional cases when really justified. The DAP measuring system is calibrated by the manufacturer and the corresponding calibration certificate is part of the delivery.

NOTICE

IMPORTANT NOTICE

THE FACTORY CALIBRATION CANNOT BE CHANGED BY USER

The user cannot change the factory calibration value of the chamber, because it is stored in the chamber memory chip. Only the DAP correction factor can be adjusted.

With the Recalibration you have the possibility to adjust the DAP for the plane in which the radiation is incident on the patient in case where absorbing materials (e.g. cosmetic cover), absorbers are permanently or temporarily present between the ionization chamber and the patient.

The procedure for determination of DAP correction factor is in general the same as which is described in Section 5.2.2 *Calibration Check Procedure*. In order to correct the deviation from Calibration Check it is necessary to read and recalculate the DAP correction factor. The new calculated factor needs to be stored by software command. For details, please refer to respective CAN EDS (Electronic Datasheet) File.

An example with a +4,51% deviation (see Chap. 5.2.2 *Calibration Check Procedure*):

Readout DAP correction factor: 1.00 (Default Value)

KermaX-plus® deviation factor: $\frac{100\%}{(100\%+4,51\%)} = 0.957$

new DAP correction factor: $0.957 \times 1.00 = 0.953$

NOTICE

IMPORTANT NOTICE

AIR DENSITY CORRECTION

When using a Ion chamber as a reference system, the reference value need an air density correction (section 6 *Air Density Correction*).

NOTICE

IMPORTANT NOTICE

DEVIATIONS TO EXPECTED MEASUREMENTS

If you observe significant deviations from the expected results please verify that no additional absorption is located in the beam (e.g. table top, additional filters). In addition please verify that the DAP correction factor of the KermaX-plus® 120-160 is set according to your system design.

⚠ CAUTION

CAUTION

DATA EVALUATION RESPONSIBILITY

The KermaX-plus® 120-160 is integrated into a system and used for data collection of Dose Area Product (DAP), the manufacturer carries the full and entire responsibility to critically evaluate the results of each measurement.

This page is intended to be empty.

6. Air Density Correction

The DAP meter is available with an built-in sensor (temperature and air pressure) for automatically correction of air density fluctuation in the ionization chamber.

If this functionality is not enabled, the user has to follow the procedure described in following Chapter.

6.1. Air Density Fluctuation in the Ionization Chamber

The ionization current is proportional to the density of air in the measuring volume of ionization chamber. According to the thermal equation of state for ideal gases, the Air density ρ is proportional to the quotient of pressure p and temperature T . It follows immediately, the correction for the influence of air density:

$$k_p = (\rho_0 / p) \cdot (T / T_0)$$

The reference values of air pressure and temperature are:

$$\begin{aligned} p_0 &= 1013\text{hPa} \\ T_0 &= 293,15\text{K} \end{aligned}$$

If necessary, using k_p to the DAP reading correction:

$$D_{dap} = M_{dap} \cdot k_p$$

For example, at a temperature of $T = 23^\circ\text{C}$ (296,15K) and an air pressure of $p = 990\text{hPa}$ is the correction of the displayed DAP M_{dap} by 3.4%.

6.2. Correction Overview

Due to built-in air-density correction sensor and variable reference sensors, air density correction for the chamber and / or the reference detector may be indicated for best accurateness.

	DAP Chamber		enabled built-in sensor	disabled built-in sensor
Reference detector				
Semiconductor detector			no action	correct DAP chamber value ¹
Ionization chamber			correct reference chamber value ¹	correct DAP chamber value ¹ correct reference chamber value ¹

¹ Temperature measurement should be done as close as possible to the corresponding chamber

This page is intended to be empty.

7. Service and Technical Support

7.1. Warranty

7.1.1. Warranty, General

IBA Dosimetry GmbH warrants the KermaX-plus® system and the associated accessories to be immaculate during the 12 months of warranty. In case of defects, IBA Dosimetry GmbH will decide whether to repair or to exchange the system.

The following information must be attached in order to claim for warranty:

- an error description
- date of purchase
- the model number of the device
- the serial number of the device

The costs for packing and shipping the goods back to the manufacturer are at the owner's expenses unless otherwise agreed upon in writing prior to shipment.

Any damages caused during the shipment of the returned goods are at the owner's expenses. On receipt of the goods the manufacturer will carry out an inspection to confirm the problem is covered by warranty conditions. If the case is not covered, the owner will receive a quote for repair.

The terms of this warranty do not affect the owner's statutory rights under applicable national or local legislation or claims against a local supplier arising from a sale or purchase contract. Claims of this nature should be addressed to the supplier.

Enhancements, modifications or repairs may only be done by IBA Dosimetry GmbH or by trained personnel with sufficient knowledge and skill who are authorized by IBA Dosimetry GmbH.

IBA Dosimetry GmbH is not liable for damages caused from the use or operation, which is not approved by the manufacturer and not described in this manual.

7.1.2. Warranty, Limitations

The warranty does not cover the following:

- Periodic checks, calibrations or preventive maintenance.
- Defects which result from modifications without the manufacturer's written approval.
- Damage resulting from normal wear.
- Damage resulting from improper use or handling including, but not limited to, the dropping or the incorrect installing of the product.
- Accidents, damages or disasters which are beyond the manufacturers control including, but not limited to lightning, fire, public disturbances and improper ventilation.
- Damages due to the transportation of the system.
- Incorrect cleaning

NOTICE

IMPORTANT NOTICE

RESPONSIBILITY, PROFESSIONAL PRACTICE

It is the OEM's responsibility to follow the guidance given in this OEM manual and to use established professional practice to ensure the product is in usable condition and is being used correctly. This manual is not intended for the use by end users.

7.2. Maintenance of KermaX-plus® 120-160

The KermaX-plus® 120-160 has no user serviceable parts which can be classified as replaceable material. The power supply and cables can be changed at any time, since this component does not have any influence on the calibration.

The KermaX-plus® 120-160 was designed to give long and reliable service and does not require special maintenance. In case the device becomes defective a repair should not be attempted but the faulty component once identified should be replaced by authorized and qualified service engineers.

CAUTION

CAUTION

REPLACING PARTS

Switch off the supply voltage and remove the external power lead from the system when replacing any parts.

NOTICE

IMPORTANT NOTICE

LOCATION OF CALIBRATION DATA

The factory calibration value is stored in the chamber memory and cannot be changed by the user.

7.3. Cleaning and Disinfection

The external parts of the KermaX-plus® 120-160 should be carefully cleaned with a soft, dry, dust free cloth. If necessary, you can clean the device using a soft, dust free cloth that has been immersed in a small quantity of pure alcohol. Under no circumstances should chemical cleaning supplies or other materials be used. When cleaning please observe the following:

- Switch OFF the system and disconnect the line voltage prior to cleaning or disinfection.
- Never spray the system with cleaning solutions. Do not allow fluids or cleaning solutions to seep into the system since it may cause damage to the equipment.
- All parts coming into contact with patients must be cleaned before each use. Clean the parts with a damp cloth, or cleaning, use pure alcohol.
- Do not use abrasive cleaners, organic solvents, and cleansers containing spot removers, etc. due to possible material incompatibility.
- Disinfectant sprays should generally not be used. The spray can seep into the system and safety features can no longer be guaranteed.
- Sprays could cause damage to electrical parts or create a flammable air / vapour mixture.

- Phenol based disinfectants and chlorine releasing preparations can weaken materials and are generally not recommended.

7.4. Troubleshooting

The following troubleshooting table may help in resolving problems.

Problem	Possible Causes
Interface is not functioning	1. Check the inter connections by disconnecting and reconnecting the cables. 2. Check the communication with another computer and Interface application.
Wrong DAP value indicated	1. Check whether the X-ray beam is penetrating the chamber correctly. 2. Check the field size setting. 3. Remove additional filtration. 4. Compare indicated values with independent dosimeter as described in the verification procedure. 5. If none of the above remedies help return the chamber to the factory for service and repair.
No DAP count during irradiation	1. Check collimator opening. 2. Check if chamber is connected and powered.
Random counting without radiation	1. Check for humidity in chamber. Dry for 24 hours and retry. 2. Check for electrical interference. If necessary call service engineer.
Test value changes	1. Check if ambient temperature is too high or too low. 2. Check once more when normal working temperature is reached.
Calibration error	1. Check for additional absorber in the beam and remove. 2. Check if the same beam quality as for the previous calibration was used.

7.5. Safety Inspection

There are no measures necessary for maintaining basic safety for the expected service life. This device does not require safety inspection.

7.6. Disposal and Recycling

This measuring device and its components contain electronic modules.

RECYCLING

None of the parts mentioned above may be disposed by the general house or hospital waste disposal system. Please observe all local regulations governing the disposal of your system. The end-user is responsible for complying with all local regulations regarding the removal of the product from service. To avoid environmental damage and/or injury and if you have no facility to convert the device to electronic waste at the end of its life cycle, please return it to IBA Dosimetry GmbH. We will ensure an environmentally correct recycling.

In the interest of complying with legal requirements concerning the environmental compatibility of our products (protection of natural resources, avoidance of waste) we endeavor to reuse components and to return them to the production cycle. We guarantee the functioning, quality and life of these components by taking extensive quality assurance measures, just as for factory-new components.

The KermaX-plus® 120-160 contains materials which can be detrimental to the environment. For this reason, the proper disposal of the KermaX-plus® 120-160 according to the regulations of national law and in view of protecting the environment must be ensured. This is guaranteed by returning the KermaX-plus® 120-160 to the manufacturer.

Materials Used:

Part	Material
Electronic	FR-4
Shielding of electronics	Zinc coated frame
KermaX chamber material	Polycarbonate; Polyphenylenether; Carbon fiber
KermaX inside frame	Aluminum, Copper
KermaX outside frame	Polycarbonate; Carbon fiber
Screws	burnished steel

7.7. Technical Support

If you need technical support, please contact the local IBA Dosimetry GmbH Service Department:

Europe, Middle East, Africa

IBA Dosimetry GmbH
Service Department
Bahnhofstrasse 5
DE-90592 Schwarzenbruck
Germany

Phone: +49 9128 607 911
E-Mail:
service-emea@iba-group.com

USA, Canada, Latin America

IBA Dosimetry America
3150 Stage Post Drive
Suite 110
Bartlett, TN 38133
USA

Phone: +1 786 288 0369
E-Mail:
service-usa@iba-group.com

APAC, Australia, New Zealand

IBA China
Dosimetry Department

No. 6, Xing Guang Er Jie,
Beijing OPTO-Mechatronics
Industrial Park (OIP)
Tongzhou District, Beijing 101111
P.R. China

Phone: +86 10 8080 9288
E-Mail:
service-apac@iba-group.com

The Quality Management system of IBA Dosimetry GmbH includes a routine to handle any reported complaints.

The user should report all complaints about the system directly to the address above.

7.8. Reporting Complaints

The Quality Management system of IBA Dosimetry GmbH includes a routine to handle any reported complaints.

The user shall report all complaints about the system to any representative of IBA Dosimetry GmbH or directly to the address above.

7.9. Returning device for repairs

Procedure for shipping the device to the factory for repair:

- Contact us via e-mail or phone at the above address to explain the problem
- We'll register your request in our system and generate an RMA form (Return Material Authorization) number
- Sending the form via Email or post would be our pleasure
- Return the parts listed on the form to us
- We will then inform you immediately that your goods have arrived and then provide an estimated repair timeframe
- A cost estimate will be generated upon request and the product will be sent back upon approval

This page is intended to be empty.

8. Technical Specifications

CAUTION

CAUTION

MAXIMUM DOSE SURFACE PRODUCT RATE

The maximum dose surface product rate of 6500 $\mu\text{Gym}^2/\text{s}$ must not be exceeded.

8.1. Specifications of KermaX-plus® 120-160 EAS

Item	Characteristic
Radiation quality	40 - 150 kVp
Energy dependence (40 – 150 kV)	$\pm 8 \%$
Quality equivalent filtration ➤ Al attenuation equivalent ➤ Al HVL equivalent (Methode of Measurement IEC 60522-1)	0,18 mm Al / 70 kV / 2 mm Al 0,10 mm Al / 70 kV / HVL 2,5 mm Al
Air Kerma rate range (at the position of the Chamber)	12 $\mu\text{Gy/s}$ - 2 Gy/s
DAP measurement range	0.1 – 42949672,94 μGym^2
DAP resolution	0.01 μGym^2
DAP linearity	$< 10 \mu\text{Gym}^2$ 5%; $\geq 10 \mu\text{Gym}^2$ 2%
DAP rate measurement range	0.10 – 6500 $\mu\text{Gym}^2/\text{s}$ max. peak [1Hz, 3.2ms radiation pulse] 13000 $\mu\text{Gym}^2/\text{s}$
DAP rate resolution	0.01 $\mu\text{Gym}^2/\text{s}$
DAP rate linearity	$< 1 \mu\text{Gym}^2/\text{s}$ 5%; $\geq 1 \mu\text{Gym}^2/\text{s}$ 2%
Min. Irritation time	$> 0.1 \text{ ms}$
Power supply	+ 15 – + 24 V _{DC} $\pm 20\%$, $< 100 \text{ mA}$
Power consumption	$< 3 \text{ W}$
Stabilization time	10 min
Chamber voltage	410 V $\pm 10\%$
Active area	Circular area with a diameter of 115 mm
Field size	69 mm ² – 9500 mm ²
Electrode spacing	6,5 mm
Communication interface	CANopen, 125/500 kBit/s

8.2. Specifications of KermaX-plus® 120-160 MFD

Item	Characteristic
Radiation quality	40 - 150 kVp
Energy dependence (40 – 150 kV)	± 8 %
Quality equivalent filtration ➤ Al attenuation equivalent ➤ Al HVL equivalent (Methode of Measurement IEC 60522-1)	0,30 mm Al / 70 kV / 2 mm Al 0,15 mm Al / 70 kV / HVL 2,5 mm Al
Air Kerma rate range (at the position of the Chamber)	12 µGy/s – up to 2 Gy/s
DAP measurement range	0.1 – 42949672,94 µGym ²
DAP resolution	0.01 µGym ²
DAP linearity	< 10 µGym ² 5%; ≥ 10 µGym ² 2%
DAP rate measurement range	0.10 – 6500 µGym ² /s max. peak [1Hz, 3.2ms radiation pulse] 13000 µGym ² /s
DAP rate resolution	0.01 µGym ² /s
DAP rate linearity	< 1 µGym ² /s 5%; ≥ 1 µGym ² /s 2%
Min. Irritation time	> 0.1 ms
Power supply	+ 15 – + 24 V _{DC} ± 20%, < 100 mA
Power consumption	< 3 W
Stabilization time	10 min
Chamber voltage	410 V ± 10%
Active area	Circular area with a diameter of 115 mm
Field size	69 mm ² – 9500 mm ²
Electrode spacing	6,5 mm
Communication interface	CANopen, 125/500 kBit/s

8.3. Specifications of KermaX-plus® 120-160 LFD

Item	Characteristic
Radiation quality	40 - 150 kVp
Energy dependence (40 – 150 kV)	± 8 %
Quality equivalent filtration ➤ Al attenuation equivalent ➤ Al HVL equivalent (Methode of Measurement IEC 60522-1)	0,30 mm Al / 70 kV / 2 mm Al 0,15 mm Al / 70 kV / HVL 2,5 mm Al
Air Kerma rate range (at the position of the Chamber)	12 µGy/s -2 Gy/s
DAP measurement range	0.1 – 42949672,94 µGym ²
DAP resolution	0.01 µGym ²
DAP linearity	< 10 µGym ² 5%; ≥ 10 µGym ² 2%
DAP rate measurement range	0.10 – 6500 µGym ² /s max. peak [1Hz, 3.2ms radiation pulse] 13000 µGym ² /s
DAP rate resolution	0.01 µGym ² /s
DAP rate linearity	< 1 µGym ² /s 5%; ≥ 1 µGym ² /s 2%
Min. Irritation time	> 0.1 ms
Power supply	+ 15 – + 24 V _{DC} ± 20%, < 100 mA
Power consumption	< 3 W
Stabilization time	10 min
Chamber voltage	410 V ± 10%
Active area	Circular area with a diameter of 144,5 mm
Field size	69 mm ² – 16300 mm ²
Electrode spacing	6,5 mm
Communication interface	CANopen, 125/500 kBit/s

8.4. Environmental Conditions and Requirements

Operational Condition	Range
Ambient temperature	+10°C – +70°C
Atmospheric pressure	700 hPa – 1060 hPa
Relative humidity range	20% – 75% (without condensation)
Transportation and Storage Conditions	Range
Ambient temperature	-20 °C – +70 °C
Atmospheric pressure	700 hPa – 1060 hPa
Relative humidity range	10% – 95% (without condensation)
Electromagnetic Conditions	Range
electromagnetic environment	professional healthcare facility environment, details chapter 2.6

8.5. Dimensions and Weight

8.5.1. KermaX-plus® 120-160 EAS / MFD

Item	Value
Dimension (L x W x H) [mm]	Ø181,3 x 18
Weight [g]	225


Figure 5 Drawing mechanical dimensions 120-160 EAS


Figure 6 Drawing position mounting holes & cable outlet 120-160 EAS

8.5.2. KermaX-plus® 120-160 LFD

Item	Value
Dimension (L x W x H) [mm]	Ø213 x 18,1
Weight [g]	279


Figure 7 Drawing mechanical dimensions 120-160 LFD


Figure 8 Drawing position mounting holes 120-160 LFD

8.6. Pin Assignment for the RJ45 8 pin Plug


Figure 9 RJ45 Plug

PIN	Signal	Description
1	CAN_H	CAN Bus Line high
2	CAN_L	CAN Bus Line low
3	CAN_GND	Ground CAN Bus / DAP Chamber power supply
4	MODE BIT 0	Node ID/Baud Rate configuration
5	MODE BIT 1	Node ID/Baud Rate configuration
6	MODE BIT 2	Node ID/Baud Rate configuration
7	MODE BIT 3	Node ID/Baud Rate configuration
8	CAN_V+	Positive power supply

8.7. CAN Node ID Configuration

Node ID	Baud rate [kBit/s]	Mode Bits *			
		Bit3	Bit2	Bit1	Bit0
0x1F	125	0	0	0	1
0x59	500	0	0	1	0
0x1B	125	0	0	1	1
0x5B	500	0	1	0	0
0x1E	125	0	1	0	1
0x1D	125	0	1	1	0
0x1C	125	0	1	1	1
0x58	500	1	0	1	0

* negative logic, 1 = CAN_GND Potential

8.8. Typical Energy Dependence


The DAP measuring device is calibrated without any absorber (see Table Reference Conditions for Calibration, Section 5.2, Calibration Check).

In the following sections the tables are for the model types with the typical results of calibration with different absorber:

8.8.1. KermaX-plus® 120-160 EAS

Tube Voltage	added Cu Filter										
	none	0,1mm	0,2mm	0,3mm	0,4mm	0,5mm	0,6mm	0,7mm	0,8mm	0,9mm	1,0mm
40 kV	1,051	0,921	0,880	0,862	0,856	0,853*	0,853*	0,853*	0,854*	0,854*	0,854*
50 kV	1,002	0,884	0,848	0,828	0,819	0,814	0,814	0,806	0,823	0,807	0,824
60 kV	0,982	0,879	0,849	0,834	0,831	0,828	0,831	0,826	0,834	0,828	0,835
70 kV	0,977	0,888	0,869	0,863	0,866	0,870	0,878	0,876	0,887	0,890	0,893
80 kV	0,981	0,907	0,899	0,902	0,912	0,921	0,934	0,937	0,950	0,957	0,964
90 kV	0,990	0,928	0,931	0,940	0,955	0,968	0,984	0,990	1,007	1,014	1,021
100 kV	1,000	0,950	0,960	0,974	0,993	1,010	1,028	1,036	1,055	1,063	1,073
110 kV	1,010	0,971	0,987	1,006	1,026	1,046	1,066	1,076	1,096	1,105	1,116
120 kV	1,022	0,990	1,013	1,034	1,057	1,079	1,099	1,111	1,131	1,141	1,152
130 kV	1,034	1,009	1,036	1,059	1,084	1,108	1,129	1,142	1,162	1,174	1,184
140 kV	1,045	1,027	1,058	1,083	1,108	1,133	1,156	1,168	1,189	1,202	1,213
150 kV	1,059	1,043	1,077	1,104	1,130	1,155	1,179	1,193	1,215	1,226	1,237


* interpolated value

8.8.2. KermaX-plus® 120-160 LFD / MFD

Tube Voltage	added Cu Filter										
	none	0,1mm	0,2mm	0,3mm	0,4mm	0,5mm	0,6mm	0,7mm	0,8mm	0,9mm	1,0mm
40 kV	0,834	0,874	0,912	0,924	0,933*	0,943*	0,952	0,959*	0,965*	0,972	0,979*
50 kV	0,894	0,937	0,963	0,970	0,972*	0,973*	0,975	0,975*	0,976*	0,976	0,976*
60 kV	0,922	0,962	0,986	0,989	0,989*	0,990*	0,990	0,990*	0,990*	0,990	0,990*
70 kV	0,941	0,978	1,000	1,005	1,004*	1,003*	1,002	1,001*	1,000*	0,999	0,999*
80 kV	0,956	0,987	1,003	1,006	1,005*	1,003*	1,002	1,000*	0,998*	0,996	0,994*
90 kV	0,961	0,990	1,003	1,003	1,000*	0,997*	0,994	0,991*	0,988*	0,985	0,982*
100 kV	0,967	0,992	1,001	1,001	0,997*	0,994*	0,990	0,987*	0,984*	0,981	0,977*
110 kV	0,972	0,992	0,999	0,997	0,993*	0,990*	0,986	0,982*	0,978*	0,974	0,970*
120 kV	0,972	0,989	0,993	0,990	0,985*	0,980*	0,975	0,971*	0,967*	0,962	0,958*
130 kV	0,971*	0,984*	0,987*	0,982*	0,977*	0,971*	0,966*	0,961*	0,956*	0,951*	0,946*
140 kV	0,967*	0,977*	0,978*	0,972*	0,966*	0,961*	0,955*	0,950*	0,944*	0,939*	0,933*
150 kV	0,960*	0,967*	0,967*	0,961*	0,955*	0,949*	0,943*	0,937*	0,931*	0,925*	0,920*

* interpolated value

This page is intended to be empty.

9. Index

§

§16 RöV	6
---------------	---

A

Absorber	27
Accessories	i, 1, 8, 17, 31
Active Area	21
Air Density Correction	29
Air Density Fluctuation	27, 29
Air Humidity	24
Air Pressure	24
Ambient temperature	40
ANSI/AAMI	3, 17
Atmospheric Pressure	40

C

Calibration	23, 46
CAN bus	1, 19
CAN Node ID	45
CISPR 11	10
Cleaning and Disinfection	32
Combination With Other Systems	8
Conventions	2
Cosmetic Cover	27
Council Directive 93/42/EEC	3
CSA	3, 17

D

DAP Correction Factor	21, 27
DAP Rate	1
Device Handling	8
Dimensions	41
Disposal and Recycling	33
Dose Area Product (DAP)	1, 5, 6, 7, 27
Dose Surface Product Rate	37

E

EC Directive 2002/96/EC	3, 18
EC Directive 2011/65/EU	18
EC Directive 93/42/EEC	6
Electromagnetic Compatibility (EMC)	9, 17
Electromagnetic Emissions	10
Electromagnetic Environment	40
Electromagnetic Environment	9
Electromagnetic Immunity	11
EN ISO 13485	17
energy dependence	46
Environmental Conditions and Requirements	40

F

Factory Calibration	26, 32
Field Size	24

H

Health and Safety Information	5
Humidity Range	40

I

IBA Dosimetry GmbH	i, 7, 8, 17, 31, 33, 34, 35
IEC 60580	17
IEC 60601	3, 6, 7, 8, 17
IEC 61000	10, 11, 12, 13, 14
Immunity to RF wireless equipment	14
Important Information	5
Installation	19
Intended Use	1
Intended User	2

L

Labels	15
--------------	----

M

Maintenance	26, 32
Materials Used	34
MDD 93/42/EEC	17
Measurement	21

O

Original Equipment Manufacturer (OEM)	i, 32
---	-------

P

Pin assignment	45
Power Supply	8, 32
Product Description	1

Q

Quality Assurance	23
Quality Management	34, 35

R

REACH	18
Recalibration	26
Recycling	33
Reference Conditions	24
Regulations, Legally Required Tests, EMC	6
Reporting Complaints	35
Responsibility	9, 19, 27
Return Material Authorization	35
RMA number	35
RoHS	18

S

Safety Inspection	33
Safety Precautions	7
Service	31
Specifications	37, 38, 39
Stability Check	23, 24
Stabilization Time	24
Symbols Used	3

T

Technical Specifications	37
Technical Support	31, 34
Temperature	24
Test tolerance and interval	26
Test Value	23
Troubleshooting	33
Tube Voltage	24

U

UL	3, 17
Unpacking	19
User Responsibility	6

V

Verification	26
--------------------	----

W


Warranty	31
WEEE	6
Weight	41

X


X-ray System	17
X-ray Tube	23
X-ray Unit	1

I hereby certify that the above copy is a true copy of the original which has been produced to me.

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