

Operational documentation for registration in Russia

Medical device: Device for hypothermia with circulating fluid TrueCool with accessories, models Premium and Lite in the composition:

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Page 1 of 68

Table of content

- 1. Description of the medical device
 - 1.1. Name of the medical device
 - 1.2. General description
 - 1.3. Intended use
 - 1.4. Indications for use
 - 1.5. Contraindications for use
 - 1.6. Professional User profile
 - 1.7. Training
 - 1.8. Improper use
 - 1.9. Medical device developer/manufacturer
 - 1.10. Risk class, classification criteria used
- 2. Using the User Manual
 - 2.1. Symbols
- 3. Safety
 - 3.1. Warnings
 - 3.2. Cautions
- 4. General system description
 - 4.1. Control unit
 - 4.2. Hoses
 - 4.3. Temperature sensors
- 5. Functional operation
 - 5.1. Operating modes
 - 5.2. Principle of operation
 - 5.3. Safety systems
 - 5.4. Water circulation system
 - 5.5. UV Disinfection system
 - 5.6. Alarm system
 - 5.7. Alarm windows
- 6. Installation
 - 6.1. Pre-installation requirements
 - 6.2. Unpacking and inspection
 - 6.2.1. Package content
 - 6.3. Moving/positioning the unit
 - 6.4. Lock/unlock the castors brakes
 - 6.5. Storage and transportation conditions
 - 6.6. Turning the device on
 - 6.7. Filling of the water reservoir
 - 6.8. Assembly of hoses and blankets
- 7. Patients preparation
 - 7.1. Sensor application
 - 7.1.1. Purpose of sensors
 - 7.2. Blankets

TrueCool

Page 2 of 68

Rev. 1.4:07

8. Operating system screens	50
8.1. TrueCool Premium – operating systems screens	50
8.1.1. Main screen description	50
8.1.2. Setting screen description	50
8.1.3. Database description	51
8.1.4. Therapy screen description	51
8.1.5. Therapy program screen	52
8.1.6. Summary screen	53
8.2. TrueColl Lite – operating system screens	53
8.2.1. Main screen description	53
8.2.2. Therapy program screen	54
9. Description of patient therapy programs	54
9.1. Manual water temperature control (Premium and Lite)	54
9.2. Fast automatic (Premium and Lite)	55
9.3. Automatic variable gradient (Premium and Lite)	56
9.3.1. Automatic variable gradient + delayed boots	57
9.4. Temperature change per hour (Premium only)	58
9.5. Therapy ending	59
9.5.1. TrueCool Lite	59
9.5.2. TrueCool Premium	59
10. Alarms and information windows	60
11. General user maintenance/cleaning principles	62
11.1. Cleaning and disinfection	62
11.2. Maintenance of the water reservoir	63
11.2.1. Draining the reservoir	63
11.2.2. Fluid circuit disinfection	63
11.2.2.1. Quarterly	63
11.2.2.2. Yearly	64
11.3. Maintenance of hoses	64
11.4. Maintenance of temperature sensors	65
11.5. Disposal	65
12. Technical specification	66
13. List of accessories	68
14. Electromagnetic interference	68
15. LED UV emitters	68
16. Shelf life	73
17. List of applicable standards	75
18. Transportation and storage of sensors and blankets	78
19. Accessories	78
Annexes	83

1. Description of the medical device

1.1. Name of the medical device

1. Control unit, in different versions Premium and Lite – 1 pcs
2. Power cable – 1 pcs
3. User manual – 1 pcs
4. Hoses - up to 100 sets (if necessary)
5. Disposable sensors, in different versions:
 - 5.1. Core temperature sensor esophageal for adults - up to 1000 pcs (if necessary)
 - 5.2. Core temperature sensor esophageal for children - up to 1000 pcs (if necessary)
 - 5.3. Skin temperature sensor for adults - up to 1000 pcs (if necessary)
 - 5.4. Skin temperature sensor for children - up to 1000 pcs (if necessary)
6. Reusable sensors, in different versions:
 - 6.1. Core temperature sensor esophageal for adults - up to 100 pcs (if necessary)
 - 6.2. Core temperature sensor esophageal for children - up to 100 pcs (if necessary)
 - 6.3. Skin temperature sensor for adults - up to 100 pcs (if necessary)
 - 6.4. Skin temperature sensor for children - up to 100 pcs (if necessary)
 - 6.5. Tympanic temperature sensor for adults - up to 100 pcs (if necessary)
 - 6.6. Tympanic temperature sensor for children - up to 100 pcs (if necessary)
7. Disposable blankets, in different versions:
 - 7.1. Blanket with clamp connections 38 cm x 53 cm - up to 1000 pcs (if necessary)
 - 7.2. Blanket with valve connections 38 cm x 53 cm - up to 1000 pcs (if necessary)
 - 7.3. Blanket with clamp connections 66 cm x 84 cm - up to 1000 pcs (if necessary)
 - 7.4. Blanket with valve connections 66 cm x 84 cm - up to 1000 pcs (if necessary)
 - 7.5. Blanket with clamp connections 66 cm x 168 cm - up to 1000 pcs (if necessary)
 - 7.6. Blanket with valve connections 66 cm x 168 cm - up to 1000 pcs (if necessary)
8. TrueCool set includes: vest, head wrap and leg pad - up to 100 sets (if necessary)
9. Head wrap - up to 1000 pcs (if necessary)
10. Vest - up to 1000 pcs (if necessary)
11. Reusable Blanket with Valve Connections 51 cm x 145 cm - up to 100 pcs (if necessary)
12. LED UV emitters - up to 10 pcs (if necessary)

Accessories:

1. Stop valve
2. Throttle
3. Mono-stereo-test-plug

TrueCool

Page 4 of 68

Rev. 1.4:07.12.2023

4. CST display 7 inch assembly
5. CST display 10 inch assembly
6. Bypass valve
7. Damper
8. Liquid level sensor
9. Pump
10. Schurter cover
11. Compressor
12. Condenser
13. Angular
14. Solenoid body
15. Solenoid coil
16. Steel frame

Target temperature management devices are used to deliberately reduce the core body temperature in patients who do not regain consciousness after return of spontaneous circulation following a cardiac arrest, as well as in newborns with hypoxic-ischemic encephalopathy. Hypothermia counteracts neuroexcitation in brain cells by stabilizing calcium and glutamate release, reducing the degree of cell death. It also stabilizes the blood-brain barrier and suppresses the inflammatory process, reducing cerebral oedema.

1.2.General description

The device consists of a system for cooling and heating as a consequences of conducted hypothermia. The medium (water) flows between the control unit and the patient's blanket(s). The essential action is achieved on the patient's body, through the direct influence of cold and heat, which affects the body's temperature. The device monitors the temperature of the human body through temperature sensors. The system actively reacts to temperature changes and adjusts the operation of the cooling-heating unit. The device is equipped with a touch-panel for operator to set work parameters.

Key features:

- Fast cooling/warming according to program and selected parameters,
- Accuracy of core temperature sensor: 0.05°C Premium and 0.1°C Lite,
- Temperature adjustment at steps of 0.05°C and 0.1°C Lite,
- Device available in two versions: TrueCool Premium and Lite,
- Easy to use 10" touch screen interface (TrueCool Premium),
- Easy to use 7" touch screen interface (TrueCool Lite),
- Core temperature (TrueCool Premium and Lite), skin temperature and tympanic temperature monitoring (TrueCool Premium),
- Multiple cooling and warming programs (automatic and manual cooling and warming programs - Premium),
- Trend screen,
- Database for off-line analysis and data export (TrueCool Premium),
- Data connection to external applications.

The medical device TrueCool achieves full operational efficiency when used with TCB/TCM blankets or a TrueCool set, manufactured by EM-MED Sp z o.o. Sp. K., which are designed for this medical device. Using any other blanket or pad on the device, could cause the intended therapy

TrueCool

Page 5 of 68

Rev. 1.4:07.12.2020

to not be carried out properly and therefore pose a risk to the patient. Furthermore it might cause damage to the device and void the warranty.

The detailed technical description of the device is included in the separate file, **Service Manual** provided by the manufacturer.

The user is obligated to comply with guidelines of relevant scientific societies in the area of therapeutic hypothermia, in terms of indications, the course of the therapy and current state of knowledge and experience of the center in a specific application

1.3. Intended use

The intended use of the TrueCool medical device is for core patient temperature management.

	Do not use TrueCool outside the facility.
	Do not use to warm/cool different equipment or food products as well as on an animals.
	Do not use the device as a place/item to sit on.

1.4. Indications for use

TrueCool is a medical device indicated for use in a healthcare facility environment by professional users only (intensive care, surgical procedures, emergency medicine).

The device can be used multiple times, with reusable blankets or single use blankets. The device is designed for single patient therapy.

Patient temperature management indications are:

- Therapeutic hypothermia,
- Inadvertent hypothermia,
- Hypoxic Ischemic Encephalopathy,
- Fever control,
- Normothermia,
- Malignant hyperthermia,
- Malignant neuroleptic syndrome,
- Stroke patients.

and other indications connected with patient temperature management. Patients:

- Infants
- Children
- Adults

Optional therapy:

The recommended therapy for Hypoxic Ischemic Encephalopathy is the cooling of the newborn after resuscitation. Childbirth-related hypoxia should begin within the first 6 hours of life: to

TrueCool

Page 6 of 68

Rev. 1.4:07.12.2020

32.5°C - 34.5°C and continue for 72 hours. After which gradual rewarming should be induced for at least 4 hours, while monitoring the child's condition.

	The warming of transdermal medications (patches) can increase drug delivery.
	During therapy pay special attention to patients with ischemic limbs and patient with arterial cross clamping as long exposure to warming may lead to burns.

1.5. Contraindications for use

Known contraindications for therapeutic hypothermia, mentioned in clinical practice guidelines, include:

severe systemic infections, confirmed multi-organ failure and pre-existing coagulopathy (fibrinolytic therapy is not a contraindication). However it is the user's responsibility to determine whether cooling/warming is appropriate for each individual patient. TrueCool is dedicated to be operated only by professional user with the medical knowledge and education, such as a medical doctor or medical personnel supervised by a medical doctor. The device should not be used on patients where clinical considerations indicate that warming or cooling of the patient is not advisable.

1.6. Professional User profile

A professional user is a registered medical facility employee (doctor or nurse), who is fully trained in the scope of conducting therapy with the use of TrueCool Temperature Management System, provided by EM-MED or its authorized representative. The user must become acquainted with this User Manual and instructions for use regarding accessories: hoses and sensors.

	TrueCool device is designed for use by, or under the supervision of, registered medical professionals only.
	User and patient cannot be the same person at the same time.
	Refer to User Manual and instructions for use regarding accessories.

1.7. Training

EM-MED or the authorized distributor of EM-MED will provide training for the system user(s) according to the intended use and operation of the device and its accessories. It is the responsibility of the hospital management to ensure that only users trained to use the equipment efficiently and safely, operate the equipment.

1.8. Improper use

Improper use of the TrueCool temperature management system can lead to skin lesions, electrical hazards, and severe changes in body temperature. The technical principles, clinical applications, and risks associated with circulatory support must be thoroughly understood before using this

TrueCool

Page 7 of 68

Rev. 1.4:07.12.2020

product. Read the entire manual before attempting to activate the system. Completion of the training program prior to using the TrueCool system is mandatory.

1.9. Medical device developer/manufacturer

Name of company: EM-MED Sp. z o.o. p. K,
Address: ul. Pocieszka 11a, 31-408 Krakow
Tel. +48 (12) 292 78 60, fax: +48 (12) 350 42 11
e-mail: office@em-med.com
internet: www.em-med.com

Technical support

Name of company: EM-MED Sp. z o.o. p. K,
Address: ul. Pocieszka 11a, 31-408 Krakow
Tel. +48 (12) 292 78 60, fax: +48 (12) 350 42 11
e-mail: techsupport@em-med.com
internet: www.em-med.com

1.10 Risk class, classification criteria used

This medical device is classified according to the nomenclature classification of the Russian Federation as Class - 2b

2. Using the User Manual

This manual describes the preparation, operation and user maintenance of the medical device TrueCool. The purpose of this User Manual is to help qualified personnel understand and operate the system and its accessories. This manual is prepared for professional personnel who may operate the system, and these operators must first have a full understanding of the proper operation of the system, thus all personnel should be familiar with all parts of this manual. If you do not understand any part of this manual, or if anything is unclear, please contact your local representative of EM-MED for further clarification. The information provided in this manual is not intended to replace regular medical training procedures. This manual should always accompany the TrueCool device. All qualified personnel operating the system should know the location of the manual and have free access, if necessity appears.

2.1 Symbols

Some symbols are used in this User Manual. The connotation of symbols is indexed below. The list encloses all markings that are visible for a user including those on TrueCool medical device, its accessories and its packaging.

	Warning: Read fully these operating instructions as they may affect safety or prevent damage		MR-unsafe! Do not expose the device to a magnetic resonance (MR) environment.
	The use of this equipment requires detailed knowledge and the observation of the Instructions for Use. This device is only intended for the described use		Do not push
	Disconnect power before maintenance activities		Indicates high temperature part or component.
	Product complies with the provisions of the Medical Device Directive 93/42/EEC as amended by 2007/47/EC.		WARNING: Make sure the UV LEDs are off before opening the device. WARNING: Use UV radiation eye and skin protection during servicing, if the UV LEDs are turned ON
	Protective earth connection (ground)		Fuse
	Type BF Equipment		USB-Connector
	Manufacturer		Refer to user manual/read the instructions
	Medical Devices		References number
	Serial number		Quantity
	Keep dry		Do not used if package is damaged
	Device is rated to a degree of Ingress Protection (IP) 2X		Keep away from sunlight

	Not for General Disposal (separate collection for electrical and Electronic equipment)		Fragile, handle with care
	Humidity limitation (10% to 70% RH)		Caution! High voltage!
	Temperature limitation (10°C-40°C)		Risk of Explosion. Do not use in the Presence of Flammable Anesthetics
	Power ON		Power OFF
	Dangerous voltage		Alternating current
	Heavy object: It is recommended that two people lift the unit to prevent injury.		Do not use handle for lift, it may cause of damage.
	Caution, electrostatic sensitive devices		For indoor use only
	Identifies the drain port		Caution
	Caution, refrigerated liquefied gases		Caution, flammable gases
	Single use		Non-sterile

	Refer to this manual for instructions and operator information.
	Read and understand all WARNINGS and CAUTIONS before you attempt to operate the system.
	Training is necessary before you attempt to operate the system.

3. Safety

3.1 Warnings

Read the entire manual before attempting to activate the system. Completion of the training program prior to using the TrueCool device is mandatory. Follow the warning notes listed below to prevent the risky situations. Follow the symbol  in the various sections of this manual to make sure that action is correct and no damage or hazardous situation occur.

	Device and all its accessories are not MRI compatible. Do not expose device to a magnetic resonance (MR) environment. The device may not function properly due to strong magnetic and radiofrequency fields generated by the MR scanner.
	Do not use the TrueCool device in the presence of flammable agents and their vapors. Risk of explosion can result – device's cooling system uses R290 refrigerant.
	The medical device contains the cooling agent R290. In case of not proper handle may lead to fire or explosion. Keep the device away from source of fire and ensured good ventilation conditions in the room. All routine maintenance should be conducted by training service technician.
	In case of using a major amount of flammable disinfecting agents, which can give off flammable vapors as well, the device should be provide with good ventilation conditions in the room.
	Do not position unit near any objects that can generate a strong electrical/magnetic field. Potential electromagnetic interference may result. Unit complies with IEC 60601-1-2.
	Do not operate the device in bathrooms or other areas where water may be splashing or may be sprayed.
	Do not allow water or other liquids to come into contact with the power connector.
	Do not use any damaged system components, accessories, cables, or sensors. This may lead to inaccurate measurements and may damage the TrueCool Temperature management System.
	Grounding reliability can only be achieved when connected to an equivalent receptacle marked "hospital use" or "hospital grade".
	Do not by-pass ground lug. Electrical Hazards may result.
	Do not connect the device's power cord to an extension cord – device must be connected directly to the facilities electrical system.
	Do not block the TrueCool device ventilation grilles. Air must be able to flow freely in and out in order to keep the device cool.
	If defibrillation is required, it is recommended to uncover the patient before the procedure.
	The vital functions of the patient should be always monitored by independent medical equipment while the therapy is carried out.

	Note that only genuine temperature sensors supplied by the manufacturer will work!
	Patient temperature should be controlled in regular intervals by the operator.
	Warming transdermal medications (patches) can increase drug delivery.
	Any liquid, substance or object present in the area between patient body and blanket/pad/mattress may influence on the therapy or patient condition due to the reaction on heat/cold conduction.
	Follow the instruction for use of the blanket/pad/mattress manufacturer in order to safety use with TrueCool device. Make sure the blanket/pad/mattress is compatible with the device and the hose nozzles fits to the connection socket of the applied part.
	The unit contain a small battery located inside the housing, therefore, the battery cannot be replaced by the user

3.2 Cautions

Follow the caution notes listed below to prevent the risky situations. Follow the symbol in the various sections of this manual to make sure that action is correct and no damage or hazardous situations occur.

	For safe handling and use of chemicals follow manufacturer guidelines.
	Users should not use cleaning or decontamination methods different from those recommended by EM-MED without first checking with EM-MED that the proposed methods will not damage the equipment.
	To ensure that the measurement system functions safely, do not subject it to high levels of humidity.
	Ensure that patient sensors/applied parts are mounted securely and that all cords and cables are well arranged to reduce the risks of injury to users, patients and damage to equipment
	Do not allow any liquids to be placed on TrueCool system (control unit / hoses or applied parts).
	Any time water is found leaking into or around the unit, connecting hose, and/or blanket, turn the unit off, disconnect the power cord from its power source, and correct the problem before proceeding.
	Proper sanitation procedures should be followed including, but not limited to, the preventative maintenance described in this manual. Leaky blankets or hoses, as well as unapproved blankets or hoses should never be used. Water leaks present a risk of infection and should be handled accordingly.
	The repair and servicing of the TrueCool system should be performed only by qualified medical equipment service technicians certified by EM-MED or authorized.

	agents of EM-MED.
	The default setting is intended to induce hypothermia/normothermia or other therapy. The system provides the physician with the option of selecting a body temperature in the range of 30°C to 40°C. If the device sounds an alarm and/or presents a display other than the standard EM-MED display, the operator should proceed according to the display message and/or the troubleshooting instructions (Chapter 10) and alarm windows (Chapter 5).
	It is strongly recommended to check the surface of the device and applied parts for freedom from mechanical damage prior to each application.
	Avoid inserting any sharp object between the patient and the blanket.
	Ensure that the initial test sound is played during start-up each time that the medical device is switched on. In case of lack of the initial test sound, user must contact with authorized service or manufacturer, as this test has a major influence on safe functioning.
	Should the therapy be stopped there will be a backfill from the pads to the reservoir, which may lead in some cases to water leakage through the overflow pipe.
	Make sure the pathways of the liquid are not blocked, bent or folded – this may affect flow between unit and applied parts.

4. General system description

One full TrueCool temperature management system consists the following parts: main unit (according to type: Premium or Lite) with power cord, 1 set of hoses (market dependent), sensors (market dependent and according to model: Premium or Lite up to 3 sensors) and optional accessories. TrueCool models differ in touch panel size: TrueCool Premium has a 10" touchscreen and the TrueCool Lite touchscreen is (7"). The TrueCool control unit is in a housing that is equipped with a handle at the front. The device is on castors so it can be moved, but only to correct the position on the working area. Two castors are equipped with brakes to prevent uncontrolled spot changing during treatment. The water fill opening is placed in the rear side (Figure 1), couplings sockets (hoses, sensors) and USB port are located on the left side (Figure 3, Figure 4). The ball valve is placed on right side (Figure 2). The ball valve is equipped with an easy to use opening mechanism to empty the tank.

	Opening the device (removal of the housing) is forbidden.
---	---

4.1 Control unit


Figure 1. Control unit description.


Figure 2. Left side.

Figure 3. Right side.

Figure 4. Rear side.


		 EM-MED Sp. z o.o. Sp. k. ul. Pocieszka 11a 31-408 Krakow, Poland			
Model	PREMIUM	Production year	2020		
Serial number	2001-00011	Voltage	220 - 240 V 50 Hz		
Weight	66 kg	Current rating	6.5 A		
Op. temp.	15 - 25°C	Input rating	1470 VA		
					
Version_#1.2_2021		Made in the European Union			
		Refrigerant R290			

Figure 5. Nameplate example(rear side).

Figure 6. Left side (Premium) Sensor's and flash drive sockets.

Figure 7. Left side (Lite) Sensor's and flash drive sockets.

4.2 Hoses

Two flexible (3,15m) hoses (1 set) are connected to the control unit and blankets to enable water to circulation. The hoses are supplied as a paired unit with one male and one female coupling to connect to the device and two female couplings to connect to blankets. Is it possibility to connect maximum 3 sets of hoses to the device, which means that maximum 3 blankets can be connected with the device. On the rear side of the device there is a retractable hook (Figure 4) which may be used for hoses (maximum 3 sets) which are not being used.

Double hose

TCH1

	Do not use the hook to place anything else on it except hoses.
	Do not bend or twist any cables or hoses. Do not use excessive force when plugging or unplugging.

Figure 8. Hoses couplings (male and female).

4.3 Temperature sensors

To monitor the patient's temperature, 3 types of sensors are used:

1. Core temperature sensor esophageal for Premium and Lite.

This is the main sensor of the device used to monitor the internal temperature of the patient's body.

2. Tympanic temperature sensor for Premium.

This is an additional sensor that should be placed in the patient's ear and is designed to monitor the temperature of the Tympanic membrane.

3. Skin temperature sensor for Premium.

This is an additional sensor that is used to monitor the patient's skin temperature.

Temperature sensors are available in two types of use: disposable and reusable, and also differ by patient type: adult and for children.

The temperature sensor is shown on Figure 10. and contains a head and a "jack" connector (plug).

Please note that only genuine temperature sensors supplied by the manufacturer will work!

Figure 9. Temperature sensor description.

Figure 10 shows a list of types of sensors and their endings (head and sockets) as well as sockets they are connected with.

Sensor type	Head	Plug (connector)	Socket
Core temperature sensor			
Skin temperature sensor			
Tympanic temperature sensor			

Figure 10. Sensor type description.

Appropriate sensor must be plugged into socket according to colours presented in Figure 10. Tympanic temperature sensor must be plugged in to the yellow tympanic socket (Premium version only), skin temperature sensor to the green socket (Premium version only) and core temperature sensor to the red socket. See Figure 10 to compare the sensors endings and the appropriate sensor sockets. Sensor sockets are located on the left side of the device (Figure 6 – TrueCool Premium, Figure 7 – TrueCool Lite).

	Do not pull cord to disconnect, always pull the connector.
	Do not immerse sensors in liquid.
	Single-use sensors are not dedicated to reuse. They have to be disposal after the use. Reuse of the single sensors may lead to faulty measurements, which in the extreme case may cause the death of the patient.
	It is recommended to connect sensors, hoses and blankets when the device is on but before a therapy program has been chosen. To read more about available therapy programs and their selection, refer to chapter 9.
	Working with electronic boards, plugs, and cables requires delicate handling.

5.Functional operation

The device consist of a system for cooling and heating the medium flowing between the central unit's water reservoir and applied parts (blankets/pads/mattresses). The essential action is achieved on the patient's body, through the direct influence of cold and heat, which affects the body's temperature. The medium circulates to and from blankets via a water pump which forces the flow.

The device provides manual and automatic programs for patient temperature management:

A) Manual water temperature control program

The software monitors the temperature of the water through (device internal) temperature sensor in the water reservoir. The software actively reacts to temperature changes and adjusts the operation of the cooling-heating unit.

B) Automatic programs

The software monitors the temperature of the human body through (double) core temperature sensor. The system actively reacts to temperature changes and adjusts the operation of the cooling-heating unit. In TrueCool Premium, the software can additionally monitor the temperature of the patient skin and/or tympanic temperature, but those temperature are not relevant to the control algorithms.

The details of each program is given in Chapter 9 of this User Manual.

5.1 Operating modes

In general, the software can indicate four operating modes of actual status of the device.

1. The STOP mode ("STOP") - no program is started or the program has been stopped.
2. The Fallback mode ("SAFEMODE"), which takes control when any of the automatic

program is turned on and no signal (or wrong reading) is received from the patient core temperature. The fallback mode operates on the principle of the manual program, where the desired temperature of the water is automatically set by the software to the value of the water temperature recorded at the moment of detection of faulty state of patient core temperature sensor. The fallback mode remains switched on till the therapy is stopped or the proper signal from the patient core temperature is received.

3. Warming mode ("WARMING") – in manual program, turned on when the desired temperature of the water is higher than measured one, at the moment of the change/set desired value of the water temperature; in automatic programs, turned on when the desired temperature of the patient is higher than measured one, at the moment of the change/set desired value of the patient temperature.
4. Cooling mode ("COOLING") – in manual program, turned on when the desired temperature of the water is lower than measured one, at the moment of the change/set desired value of the water temperature; in automatic programs, turned on when the desired temperature of the patient is lower than measured one, at the moment of the change/set desired value of the patient temperature.

Operating mode indication is shown next to the "STATUS" text on the therapy screen. Additionally, the Premium version indicates the mode by changing the background color on the therapy screen (blue – cooling, orange – warming).

The device cannot change the operation mode by itself, the only exception is to going into the SAFEMODE from

"WARMING" or "COOLING" mode in automatic programs.

The implemented alarm system constantly monitors the actual mode of the device and in case of any deviation from the normal or desired state informs the user by appropriate visual/audio indication.

Temperature values averaging.

Any patient temperature value shown on the device's screen is not subjected to any averaging - real temperatures from the sensors are shown on the screen with no delay coming from the software.

5.2 Principle of operation

Manual water temperature control program

The program's operation principles are based on controlling the water temperature in the water reservoir. The user defines the target water temperature and the appropriate functions of the components are used to achieve and maintain desired temperature within stated range. Those functions are warming (thanks to the heater installed in the bottom water reservoir) and cooling (thanks to the compressor, electromagnetic valve which controls the compressor and copper pipes present in the bottom water reservoir to provide heat absorption from the water, the heat is distracted due to the condenser with rotating fan).

Automatic programs

The program's operation principles are based on achieving the desired patient's core temperature by controlling the appropriate functions of the components (warming and cooling) in a way so the defined conditions of the algorithms are maintained (see Chapter 9). Automatic programs use the assumptions of the Manual water temperature control program, but the water temperature value is set automatically by the software, based on the particular algorithm assumptions, so the change in patient core temperature is achieved (and maintained). Algorithms of the automatic programs need

a feedback variable from a patient which represents the condition of the patient or the influence/efficiency of the applied therapy on the patient.

The feedback variable used in the implemented algorithms is patient core temperature calculated by measuring the rate of resistance changed due to the patient core temperature change. The expected latency/delay times are within hundreds of nanoseconds with no noticeable delay for the operator. It is very important to provide accurate readings of the patient core temperature, so the algorithm can apply appropriate functions, this is why the placement of the temperature sensor is crucial and any disturbance to the core temperature must be minimized.

Disturbances may affect the result coming from the energy exchange between patient and the blanket(s). Sources of disturbance may come from the patient (change in temperature due to infusions) or from outside the patient (environmental conditions such as room temperature and humidity). The influence of this disturbance on the patient temperature is a function of several factors (like weight of a patient, effectiveness / quality of how well the patient is wrapped in blanket/pad, amount and temperature of infusion liquids or other drugs, which also may be different from patient to patient. These factors or combinations thereof may also vary with time, therefore we end up with numerous possible influences of varying importance in which can impact the patient temperature.

It can be assumed that the highest possible importance of the impact may occur in case of small patients, where the difference between desired and actual patient temperature is relatively small. In that case it may happen that the disturbance can "easier" "push out" the values of the patient temperature from the programmed that is to be maintained. This is limited if the disturbance is constantly present (like constant room temperature), but can be significant in case of rapid changes (like infusion of a pre-warmed solution). On the other hand, with an adult patient fully covered by the blankets/pads, through which cold water circulates cooling the patient, the influence of the constant room temperature is minor, and also infusion of pre-warmed liquids may have lower impact on the patient temperature than in the case of neonates. Of course, it also depends on the amount of the liquid to be administered so the comparison of the impact on different patients may not be accurate. The influence of the disturbances, especially rapid ones may have a positive impact on patient temperature (or rather the change of the patient temperature) when used on purpose and thoughtfully – at the stage of the therapy where fast cooling is necessary, infusions of liquids (not necessary pre-cooled, but at least at room temperature) may speed up the change in the patient temperature, but the same infusion may be considered as unwanted in case, when we need to increase the patient temperature to a normal level. The same disturbance, depending on the operation mode can be useful or cause a danger to the patient.

Due to the human variation resulting from the wide range of the patients (neonates, children and adults) and also different influence of potential disturbances on the patient, an appropriate automatic therapy program should be used.

Automatic programs expect the feedback variable in a specific ranges. Range of 0-50°C assumes a proper reading from the sensor, while a range of 28-41°C correspond to a reasonable values of the patient temperature. The alarm system defines the appropriate alarm condition to inform the user if the value is outside defined ranges and switches to the SAFEMODE if necessary. Additionally, patient core temperature sensor consists of two independent resistance elements measuring the change in resistance – if one value is significantly different than the other one, the alarm system will indicate the possible damage of the core sensor.

The patient's core temperature (feedback variable) is the only source needed for automatic operation. It's value, depending on the chosen program, influences the rate of energy which is to be delivered or absorbed to or from the patient. This transfer of energy is in the form of heated/cooled water circulating in the applied part (pad/blanket) which is then transferred to the patient by direct contact (conduction).

The main goal of the device is to reduce, maintain or slowly increase the patient's body temperature which is done through the conduction of applied parts (blankets/pads) and patient body, where the applied parts have a specific temperature. The change in patient's body temperature (feedback variable) depends on:

- temperature difference between surface of the pads and patient temperature,
- the surface area of contact between pads and patient body,
- the weight of the patient (new-born, child or adult),
- environmental conditions (room temperature, humidity),
- other factors like infusions of pre-warmed or cool liquids).

By applying the same therapy conditions to two extremely different patients, a different response can be observed and therefore this type of device must be constantly supervised by a professional user.

Because of these difficulties appropriate algorithms (programs) were implemented in the software so that the user, by choosing the appropriate parameters of a particular program can control the response of the system and minimize unwanted influence of factors affecting the therapy.

General assumptions to the alarm system were implemented and can be referred to the relevant parameters of the system response idea:

Settling and response time

Response time is time when patient's temperature reaches 90% of the set temperature and patient's temperature is still moving to the set temperature. The settling time is a period when patient's temperature reaches 95% of the set temperature and the difference between patient's temperature and the set temperature does not exceed $\pm 5\%$ of value of set temperature.

In general the importance of settling and response time can be split into three situations. First, when we need cooling as fast as possible, second, when we need to slowly increase the temperature as too fast warming may lead to death of a patient, third, when we need to adjust the temperature of a patient.

In case of too fast rewarming of the patient we use Medium Priority Alarm no. 1 - Dangerous increase of the patient temperature.

For other cases we use Low Priority Alarm no. 1 - No cooling or warming of the patient,

which informs the user that the patient temperature has not changed whilst the change was expected.

The settling and response time depend on the chosen algorithm and its parameters.

Relative overshoot, Command overshoot, Steady state deviation

Command overshoot is the maximum positive deviation when patient's temperature was out of acceptable limit of +5% of value of the set temperature.

Relative overshoot is the difference between the maximum positive deviation when patient's temperature was out of acceptable limit of +5% of value of the set temperature and the final patient's temperature, which is in the acceptable limit of +/-5% of value of the set temperature.

Steady-state deviation is the difference between the set temperature and the lower of the upper limit of acceptable difference between the patient's temperature and the set temperature. An overshoot may not be observed in the case of adults, but may be apparent in the case of neonates. This is due the "inertia" of the patient's body, which is relatively "small" in case of adults, and "high" in case of low-weight patients. In the other words small patients may quickly loose body temperature especially when the heat exchange between blanket and body is efficient (high temperature difference). In those cases, two alarms inform the operator about possible dangerous situations: Low Priority Alarm no. 2 - Patient core temperature out of the programmed range and Medium Priority Alarm no. 2 - dangerous temperature of the patient. The last one is also useful to notify the user about variations from the reference variable.

The overshoot and deviation from the reference variable depend on the chosen algorithm and its parameters.

Tracking error

Tracking error is the difference between patient's temperature and set temperature at the same time during the therapy.

The tracking error will vary depending on the patient's body and effectiveness of the heat exchange. The device is capable to achieve maximum and minimum value of the manipulated variable within 20 minutes (but this also depends on many factors like the number of blankets attached, environmental conditions) and in 20 minutes the reference variable can change by a different value depending on the patient and therapy conditions. In that time for small patients the error may be quickly reduced whilst in case of adult patient the change may be very slow.

The error depends on the chosen algorithm and its parameters.

Variable logging

All the therapy data (including patient's temperatures and TSET) is saved into storage memory and cannot be deleted (power losses does not delete the saved data). The data can be accessed by the authorized service or manufacturer. TrueCool Lite do not provide access to this data for the user, while TrueCool Premium plots the therapy data (temperatures) in real time and allows therapy data to be exported to an external storage via USB port. If the storage memory is full, the oldest data is deleted.

5.3 Safety systems

The safety systems in TrueCool devices consist of an electronic alarm system and mechanical overheating protection. The electronic alarm system is embedded in the software and constantly analyses the operation of the device and its components.

The alarm system appearance is described in details in Chapter 5.5 and protection values from overheating and overcooling are given in Chapter 12.

Independent, mechanical overheating protection is supplied by a bimetallic thermostat attached to the water reservoir which will turn off the heater if the measured temperatures reach a level 65°C.

Undercooling below 20°C is not possible due to the physical aspects of heat exchange. The following description supplements chapter 5.5 in relation to more technical aspects of defining alarm conditions.

Therapy not efficient?

Alarm is valid only for automatic programs.

Alarm is active for the first 30 minutes from the confirmation of the chosen program or modification of the program setting.

After 30 minutes, the alarm is shown or

turned off. The alarm is shown if the

difference between TPATIENT (core)

measured with an interval of 30 minutes is

less

or equal 0.1°C (condition for program no 2 and 3), and

Possibly no efficient cooling or warming of the patient. The alarm is desired to point that the patient temperature has not been changed – due to the chosen

therapy program we expect that the temperature of the patient is changed. The alarm is of the latching type

– if shown, it will not disappear without user interaction.

additionally TGRADIENT must be higher or equal 0.25°C in case of program no 4 (TrueCool Premium only).

Patient core temperature out of the programmed range.

Alarm is valid only for automatic programs.
Alarm is active when the TPATIENT (core) falls in the appropriate range:

±0.05°C from TSET for TrueCool Premium

±0.1°C from TSET for TrueCool Lite

Alarm is deactivated if TSET has been changed and the TPATIENT is out of the defined range.

Patient core temperature out of the programmed range.

If the TPATIENT value falls in the appropriate range the alarm conditions are activated. The alarm condition will be deactivated if the TSET is changed and TPATIENT value is out of the appropriate range.

If TPATIENT is out of the programmed range (as specified in the program algorithm), the alarm condition is true, until TPATIENT value will fall back to the programmed range.

The alarm is of the latching type – if shown, it will not disappear without user interaction.

Water temperature out of the programmed range.

Alarm is valid only for manual program (no 1).

Alarm is active when the water temperature falls in the range of ±0.5°C from TSET.

Alarm is deactivated if TSET has been changed and the water temperature is out of the defined range.

Water temperature out of the programmed range.

If the water temperature value falls in the range of TSET±0.5 degrees of Celsius, the alarm condition is activated. The alarm condition will be deactivated if the TSET is changed and water temperature value is out of the range of TSET±0.5 degrees of Celsius.

The alarm is of the non-latching type – if shown, it will automatically disappear if the condition of the alarm is false.

Low flow.

Conditions are active only if the program has been chosen and accepted.

Alarm conditions are turned off for 20 seconds after the therapy has been started and remains active until therapy is stopped.

Water flow is constantly measured by an internal flowmeter. If the water flow rate is less or equal than 0.5 liter per min, than the alarm condition is true.

Low flow alarm will indicate that the flow through the water circulation system or blanket itself has been compromised.

The alarm is of the non-latching type – if shown, it will automatically disappear if the condition of the alarm is false.

Replace UV LED soon!

Condition is always active. The device constantly monitors the work and time of working of UVC modules. If the counter reaches the value of 700 hours the appropriate message will be shown.

Replace UV LED soon as the UV LED efficiency is limited due to its lifetime.

The UVC time counter may be reset to zero only by a service technician after the UVC modules are replaced by new one.

The alarm is of the latching type – if shown, it will not disappear without user interaction.

Dangerous increase of the patient core temperature!

Condition is always active, if the patient core temperature sensor is connected to the device.

The alarm is shown if the difference between PATIENT measured with an interval of 1 minute is higher or equal 2°C.

Dangerous increase of the patient temperature. The condition is checked every 60 seconds (1 minute interval).

The alarm is of the latching type – if shown, it will not disappear without user interaction.

Dangerous core temperature of the patient!

Condition is always active, if the patient temperature sensor is connected to the device.

The alarm is shown if measured patient core temperature is higher than 41°C or is lower or equal to 28.0°C.

Dangerous temperature of the patient. Please check if the patient core temperature sensor is placed correctly. The alarm is of the latching type – if shown, it will not disappear without user interaction.

Patient Core Temperature Sensor Fault!

Condition is always active, if the patient temperature sensor is connected to the device.

The alarm is shown if measured patient core temperature is higher than 50°C, is lower or equal to 0°C or no signal from the sensor is received.

Fault of the patient core temperature sensor.

If the core temperature sensor is connected and is faulty, then the error will be displayed (also, instead of temperature value the "ERROR" text will be shown).

If the core temperature sensor is not faulty, but it has been disconnected during the therapy (excluding manual water temperature control program), the error will be displayed and the device will change the operation mode to fallback mode (indicated as the SAFEMODE).

The alarm is of the latching type – if shown, it will not disappear without user interaction.

5.4. Water circulation system

Water is pumped in and out of the blanket using a water pump.

5.5. UV Disinfection system

A UV disinfection system is installed in the dual water reservoir which works in an automatic 4-hour a day cycle irrespective of a device to patient connection, hindering pathogen growth. The UV disinfection system has been tested and its effectiveness has been confirmed by tests performed by an independent laboratory.

5.6. Alarm system

In the device an alarm system conforming to the standard EN 60601-1-8 is used. In the following section the specifications and basic settings are described.

Alarm classification

Alarms are classified regarding the type and risk as:

- Physiological Low Priority Alarms.
- Physiological Medium Priority Alarms.
- Technical Low Priority Alarms.
- Technical Medium Priority Alarms.

Alarms with high priority are not used, because the device belongs to class IIb and, therefore it is not intended for the monitoring of vital physiological parameters, where the nature of variations is such that it could result in immediate danger to the patient.

The software generates distinctive visual and auditory alarm signals for each alarm priority.

Physiological alarms are activated when the software detects a potential hazardous situation based on the readings from the patient core temperature sensor. Physiological medium priority alarms must always be confirmed by the operator and even if the condition of the alarm is not valid anymore, the alarm window will remain present on the screen until the appropriate action button is pressed.

Technical alarms are activated when the software detects a probable malfunction of the instrumentation or of a connected accessory. In contrast to physiological alarms, technical alarm windows disappears if the alarm condition is no longer true, and the alarm cannot be paused or switched off until the condition is no longer true.

Visual alarm signals

Visual alarm signals appear on the left top part of the display when an alarm is active (the condition of the alarm is valid). The following visual alarm signals are used:

Alarm priority	Color	Flashing frequency	Duty cycle
I - low priority	yellow	constant (on)	100% on
II - medium priority	orange	1 Hz	50 % on

Table 1. Alarms priority description.

When a low priority alarm condition is detected, an appropriate alarm symbol is displayed in frame on the top part of the screen (excluding LPA_1 alarm). As the exception to the rule of low priority visual alarms, LPA 1, LPA 15 and LPA 16 are indicated by the window only, and this window is displayed only once (no additional indication after the alarm is off). If more information about an alarm is needed, press the alarm symbol displayed on the screen.

When a medium priority alarm condition is detected, an appropriate alarm window is displayed on the centre of the screen. The alarm window provides available action buttons with appropriate function. When the alarm window is closed (which is possible only if the operator presses the appropriate button), an additional indication of the alarm is displayed in the frame on the top part of the screen with the exception of the "critical water level" and "no flow" alarms, which are indicated on the screen by replacing the original pictograms of water level and flow.

Below each symbol displayed on the screen, text will indicate if the alarm is associated with physiological ("PATIENT") or technical ("DEVICE") alarm condition.

Auditory alarm signals

Auditory alarm signals sounds once an alarm condition is detected and every 10 seconds if the medium priority alarm was not Paused or OFF (visual alarm signal – window – is displayed on the screen).

The following melodies are used for auditory alarm signals (according to EN 60601-1-8 Annex F):

Alarm priority	Type	Notes
I - low priority	All	e c
II - medium priority	General	c c c
II - medium priority	Temperature / energy	c d e
II - medium priority	Device fault	C c c

Table 2. Audio alarms signals description.

Reminder signal

Single "beep" sound used to:

- draw the user's attention when the action is required or advised,
- indicate that the low or medium alarms were PAUSED.

Power loss

If the therapy is running, a break in power supply will cause the alarm system to turn on the audible constant signal, which lasts up to one minute.

Alarm condition and generation delay

If an alarm condition is detected, the appropriate indication of the visual and auditory alarm is applied immediately. The time from the triggering event to the generation of the alarm depends on the speed of signal transmission between components and takes less than a second. No software delay in alarm generation was implemented, except "no flow" (MPA 8) and "low flow" (LPA 6) alarms, which conditions are turned off for 20 seconds from the start of the therapy. This delay is included due to the delay in starting the pump and the time needed to force the water to circulate through the pipes and flowmeter.

5.7. Alarm windows

Structure of the alarm window

1. Visual alarm symbol area.
2. Alarm name and description area.
3. Additional information area.
4. Buttons area.

Example of the alarm window:

Table 3. Alarm window structure.

In case of medium priority alarms, the symbol in the field marked as 1 is “blinking” as it is defined in the Table 1.

Buttons in windows

The alarm windows provides buttons, with appropriate options for the operator:

- “PAUSE” – pauses the alarm;
- “OFF” – switches off the alarm;
- “OK” – confirmation.

Some windows have an informative function and provide only the “OK” button, which is used for confirmation only

(only one option is available – no choice for the operator).

In some cases, the window will also provide a “CANCEL” button, which is used to close the window – no action or operator choice is required.

In one case a special type of button is used for suspending the alarm system for physiological conditions (in the reset alarm window). Buttons are shown in the table below.

Button name	Symbol	Description
OFF		Visual alarm signal behaviour after the button is pressed 1. The alarm window is closed. 2. If the condition of the alarm is true, the appropriate on-screen visual alarm signal is displayed. 3. The Alarm window will not be displayed again.


		Auditory alarm signal behaviour after the button is pressed 1. If the condition of the alarm is true, the Single Auditory Signal sounds every 60 seconds.
PAUSE		Visual alarm signal behaviour after the button is pressed 1. The alarm window is closed. 2. If the condition of the alarm is true, the appropriate on-screen visual alarm signal is displayed. 3. After 5 minutes (from pressing the action button) the alarm window is displayed again. Auditory alarm signal behaviour after the button is pressed 1. If the condition of the alarm is true, the Single Auditory Signal sounds every 60 seconds. 2. After 5 minutes (from pressing the action button) the Low Priority Auditory Signal sounds.
OK		The (alarm) window is closed. If the action/function was implemented for the button, the action is triggered. If no action/function was implemented, no action is triggered.
CANCEL		The (alarm) window is closed. No additional action is triggered.
SUSPEND		The (alarm) window is closed. The appropriate visual indication is displayed on the screen.

Table 4. Buttons description.

Priority	Low	Medium
Alarm Condition		
Alarm OFF		
Alarm OFF Acknowledged		
Alarm Paused		
Alarm Paused Acknowledged		
Audio Alarm Paused Acknowledged	 or 	

Audio Alarm OFF Acknowledged	 or 
Audio Alarm OFF Used as a button to turn audio alarm off.	
Audio Alarm Paused Used as a button to pause audio alarm.	

Table 5. Symbol description.

Below each symbol text will indicate if the alarm is associated with a physiological (PATIENT) or technical (DEVICE) condition.

Screen indications of alarms – resulting from the chosen action button from alarm windows

The screen indication alarms are displayed on the top part of the screen in a frame. The frame with list of the alarms (icons for each detected condition) is shown when the alarm icon is pressed.

Action button	Screen indication for the physiological alarms		Screen indication for the technical alarms	
	Premium	Lite	Premium	Lite
PAUSE				
OFF				
Action button	Screen indication for water level		Screen indication for water flow	
	Premium	Lite	Premium	Lite
PAUSE				
OFF				

Table 6. Screen indications.

Screen indications of alarms –defined by default to be indicated on the screen.

Action button	Screen indication		Additional description
	Premium	Lite	
OK			If the on-screen visual alarm indication is pressed, then the appropriate information window from table 8 is displayed.

OK			If the on-screen visual alarm indication is pressed, then the appropriate information window from table 8 is displayed.
OK	Water level 	Water level 	Normal level of water – not full, but above medium.
OK	Water level 	Water level 	If the on-screen visual alarm indication is pressed, then the appropriate information window from table 8 is displayed.
PAUSE / OFF	Water level 	Water level 	If the on-screen visual alarm indication is pressed, then the appropriate alarm window from table 8 is displayed.
	Water flow 	Water flow 	If the on-screen visual alarm indication is pressed, then the appropriate information window from table 8 is displayed.

Table 7. On-screen visual signals description.

Configuration of the alarm windows

If the star: "*" is used at the end of the alarm name, the alarm is to be implemented only in Premium model.

Symbol	Alarm name and description.	Additional information	Buttons
	Therapy not efficient! The core temperature of the patient measured 30 minutes ago has not been changed.	Recommendation to check if: - the patient is covered properly - there is a proper flow through blanket(s) - the water temperature is correct - the temperature sensor is placed correctly.	
	Patient core temperature out of the programmed range.	This might be caused by different factors. Observe the condition of the patient. Check if no other alarm indication is presented on the screen.	OK
	Water temperature out of the programmed range.	This might be caused by different factors. Observe the condition of the patient. Check if no other alarm indication is presented on the screen.	OK
	Water level at medium level.	Add water until you hear a single beep.	OK
	Low flow.	Check if: - blanket(s) are properly connected to the hoses - hoses or blankets are not kinked - hoses are not crushed	OK
	Fault of the condenser temperature sensor.	Call service.	OK

	Fault of the secondary temperature sensor in the block.	Call service.	OK
	Fault of the patient skin temperature sensor.	Check the patient skin temperature sensor, replace sensor or call service.	OK
	Fault of the patient tympanic temperature sensor. Error Code: LPA10	Check the patient tympanic temperature sensor, replace sensor or call service.	OK
	System Error. UV LED regulation system fault. Error Code: LPA11	Call service.	OK
	Potential fault early warning system. Error Code: LPA 12	Increased risk of pump failure. Call service.	OK
	Potential fault early warning system. Error Code: LPA13	Increased risk of heater fault. Call service.	OK
	Potential fault early warning system. Error Code: LPA14	Increased risk of compressor fault. Call service.	OK
	Service reminder!	Maintenance required due to device running time. Call service.	OK
	Replace UV LED soon! Error Code: LPA 16	The UV LED diode is nearing to the end of its lifetime. Call service.	OK
	Dangerous increase of patient core temperature! Error Code: MPA1	The patient core temperature increased by 2°C in 5 minutes! Check if: - temperature sensor is placed correctly and - the condition of the patient!	
	Dangerous core temperature of the patient! Error Code: MPA2	The patient core temperature has reached an extreme value! Check if: - temperature sensor is placed correctly and - condition of the patient!	
	Primary Overheating Alarm! Water temperature above the limit! Error Code: MPA3	The heating function is turned off until the temperature of the water drops below the limit. Check patient condition. If the error message not disappear, please contact the service technician.	
	Secondary Overheating Alarm! Water temperature above the limit! Error Code: MPA4	The heating function is turned off until the temperature of the water drops below the limit. Check patient condition.	

		If the error message not disappear, please contact the service technician.	
	Primary Overcooling Alarm! Water temperature below the limit! Error Code: MPA5	The cooling function is turned off until the temperature of the water rises above the limit. Check patient condition. If the error message will not disappear, please contact the service technician.	
	Secondary Overcooling Alarm! Water temperature below the limit! Error Code: MPA6	The cooling function is turned off until the temperature of the water rises above the limit. Check patient condition. If the error message will not disappear, please contact the service technician.	
	Critical Water Level! ADD WATER! Error Code: MPA7	Operation of the device may be limited due to the safety of the unit components. Check for any water leakage. Add water until audio indication is heard.	
	No Flow! Error Code: MPA8	Check if: - blanket(s) are properly connected to the hoses - hoses or blankets are not kinked, or - hoses are not crushed	
	Primary Temperature Sensor Fault! Error Code: MPA9	Call service.	
	Patient Core Temperature Sensor Fault! Error Code: MPA10	Call service.	
	Pump fault! Error Code: MPA11	Call service.	
	Heater fault! Error Code: MPA12	Call service.	
	Compressor fault! Error Code: MPA13	Call service.	
	Cooling system efficiency problem.	Call service.	
	High condenser temperature!	Call service.	


	Fault of the mechanical water level sensor!	Call service.	
	System Error! Error Code: _____	Call service.	



Table 8. Alarm windows description.

Reset of alarms – window

Reset alarm window will appear on the screen, if one of the on-screen visual alarm signal indication is pressed for at least 3 second.

Alarm Reset Window:

Action button	Screen indication		Description
	Premium	Lite	
 RESET CHOSEN ALARM	N/A	N/A	If the button is pressed: the chosen alarm is turned off, the alarm reset window is closed, the alarm system will re-check the condition(s) for the chosen alarm and will trigger the appropriate action (alarm) if the alarm condition(s) is true. If audio was OFF or Paused (indicated by the on-screen visual signal), then the Audio Pause or OFF function is cancelled. If alarms were suspended (indicated by the on-screen visual signal), then suspension of the alarms is cancelled.
 RESET ALL ALARMS	N/A	N/A	If the button is pressed: all the alarms are turned off, the alarm reset window is closed, the alarm system will re-check the condition(s) for all the alarms and will trigger the appropriate action (alarm) if the alarm condition(s) is true. If audio was OFF or Paused (indicated by the on-screen visual signal), then the Audio Pause or OFF function is cancelled. If alarms were suspended (indicated by the on-screen visual signal), then suspension of the alarms is cancelled.


 AUDIO PAUSE (ALL ALARMS)			If the button is pressed: all the auditory alarm signals are turned off for 10 minutes, the alarm reset window is closed, the appropriate on-screen visual alarm signal is displayed, after 10 minutes all the auditory alarm signals are turned on, the appropriate on-screen visual alarm signal is turned off. If the on-screen visual alarm indication is pressed, then the Alarm Reset Window is displayed.
 AUDIO OFF (ALL ALARMS)			If the button is pressed: all the auditory alarm signals are turned off permanently, the alarm reset window is closed, the appropriate on-screen visual alarm signal is displayed. If the on-screen visual alarm indication is pressed, then the Alarm Reset Window is displayed.
 SUSPEND ALARMS FOR 10 MINUTES	 PHYSIOLOGICAL ALARMS SUSPENDED!		If the button is pressed: all the alarm signals are turned off for 10 minutes, the alarm reset window is closed, the appropriate on-screen visual alarm signal is displayed. If the on-screen visual alarm indication is pressed, then the Alarm Reset Window is displayed. If the button is pressed when the alarms has been already suspended, the alarm suspension time is counted from the beginning.
 CANCEL BUTTON	 N/A	 N/A	If the button is pressed the reset alarm window is closed with no additional action and indication.

Table 9. Alarm reset window description.

Auditory alarm indication rules

The auditory alarm signals are generated at the same time when the alarm window (visual alarm signal in the form of window) is displayed.

Overlapping of multiple auditory alarm conditions

In case when the multiple conditions of the different auditory alarm indications are met, the following rules applies:

- different priorities of the auditory alarm: the highest priority alarm is generated (the "melody" of the lower priority alarm is turned off and immediately the "melody" of the higher alarm priority is generated),
- if different priority conditions are met, and the condition of the highest priority alarm changes to false, then the rules of the lower priority alarms applies (as if they were detected for the first time),
- if multiple alarm conditions of the same priority are met, in case when condition of one of the alarm is false but the other one is still met, the rules remain unchanged for generation of the given alarm priority,
- if during generation of the auditory signal, another condition of the same type and priority is detected, the generation of the alarm is not disturbed by the new alarm generation,
- if during generation of the auditory signal for an alarm of the physiological type, another

condition of the same priority is detected, the generation of the alarm is not disturbed by the new alarm generation, if the detected condition is also associated with the physiological type,

- if during generation of the auditory signal for an alarm of the other than physiological type, the other condition of the same priority is detected, the generation of the alarm is changed by the new alarm generation, if the detected condition is of the physiological type,
- in case the alarms (or auditory alarms) are paused and multiple alarm conditions are detected, the above rules are applied.

Visual alarm indication rules

The visual alarm signals are generated at the same time when the auditory alarm signal is generated. Visual alarms are displayed on top of each other. In the instance of multiple alarms the following rule applies.

Overlapping of multiple visual alarm conditions

In case when the multiple conditions of the different visual alarm indications are met, the following rules applies:

- if the multiple alarms are of the same type and priority, the windows are displayed by time of occurrence (the last window is "on top of the other(s)"),
- if the multiple alarms are of the same priority but different type, and no physiological type is detected, the windows are displayed by time of occurrence (the last window is "on top of the other(s)"),
- if the multiple alarms are of the same priority but different type, and any of the alarm is of physiological type, the window which is displayed is to be for physiological type of the alarm ("on top of the other(s)"),
- the alarm windows are displayed "above" (on top of) the other screens and (informative) windows and are always visible prior to other types of screens.

In case of multiple alarm windows which are displayed (but not visible) on the screen, if the button on the last (visible) window is pressed, the following action is performed:

- all the alarm windows for a technical conditions are closed and if it is possible the same closing function ("OFF" or PAUSE") is applied to not visible windows, otherwise the "OK" closing function is applied (as if the operator has pressed the OK button) or any other available button.
- the above functionality is not applied to physiological alarm windows, which mean that every alarm associated with the patient has to be confirmed by the user. In the case of multiple alarm windows for physiological conditions, when the one window is closed, the other one is shown (applying the rules described in the "Overlapping of multiple visual alarm conditions" section).

6. Installation

6.1 Pre-installation requirements

TrueCool has a space requirement of: 70cm x 60cm x 110cm (w x d x h). The device is dedicated for indoor use only with environment and electrical requirements as listed below:

- The room lighting must comply with the applicable health and safety regulations,
- The room temperature should be between 15°C and 25°C,
- Relative humidity in the room is 30-85%,
- The product does not require a fixed ambient atmosphere pressure,
- The control unit must be connected to electricity supply according to requirements: 230VAC (market dependent).

For Indoor use only.

6.2 Unpacking and inspection

	The unit should be unpacked, installed and tested only by authorized/trained personnel.
	The package is heavy. Two people should lift the package to prevent injury.
	Wait 24 hours before first use, if transportation to facility was conducted in low temperatures or the transportation conditions were not certain.
	If you suspect that the device was transported in a horizontal position, contact the manufacturer.
	Do not use if package is damaged. Report to the transportation company.
	Do not push.

The delivery must be unpacked according to information contained in this short instruction. During the unpacking process, look carefully for signs of shipping damage. If any unacceptable damage is found, notify the transportation company immediately and file a claim. The transportation company is responsible for the shipment after it leaves the factory. If problems other than shipping damage are found, notify your supervisor and/or EM-MED representative. Information about packaging to be kept as it may be used in the future. Check carefully, if the control unit and its accessories are not damaged in any way. TrueCool must be transported to facility in vertical position only.

Step	Description
1	Remove top foams from the device.
2	Lift cardboard and take off it from the devices. The original package has to be kept.

	
3	Two people have to catch the device by holder and the bottom of the device. Next the device should be carefully tilted approximately 20° and lifted in order to remove the bottom foam and cardboard by third person. CAUTION: It is forbidden to lift the device with holder by one person, because it may be a cause of its damage (see pendant on the device). 
4	The device has to be again placed on a flat surface and with blocked wheels to avoid uncontrolled movements of the device. Next, the device should be carefully checked, that it was not damaged during transport. If damage of the device is found, please follow the instruction described above.
5	If you have information, that transport was being conducted in low temperature, please wait 24 hours before first use.

6.2.1 Package content

Package content differ and depends on list of accessories ordered by user and TrueCool model. The typical package consist of:

- User manual,
- Device TrueCool Lite or Premium with power cord,
- Sensors (according to model 1 to 3 pieces),
- Hoses (1 to 3 sets, according to order).

TrueCool

Remove the TrueCool from service if the casing or membrane control panel is cracked or internal components are exposed. Contact with internal components could result in electric shock or thermal injury to the patient or operator.

6.3 Moving/positioning the unit

The device must be transported in line with internal facility procedures. The castors of the device are dedicated only for in room positioning. The device has a handle (see Figure 1), which is to be used for positioning of the device. The positioning procedure is conducted while the brakes are unlocked. Positioning procedure means moving forward/backward, right/left. Do not position TrueCool so that it is difficult to disconnect the device from the mains electricity supply. Do not position device in place, where air flow from control unit will be blocked, this may cause over temperature and damage of the device. Also do not position the device near any objects that can generate a strong electrical/magnetic field as this may result in potential electromagnetic interference. Position the unit so that the touchscreen faces the operator. The control unit may be positioned on flat and horizontal surfaces only. It is recommended to set the position of the control unit: screen fronted to the user or open space and hoses connection on the bed side.

If it is necessary to move the device between rooms, should be used internal transport system. The device is not intended for passing through thresholds and jumps. In case not to follow the above advise, it might lead to damage of castors and/or the device.

The maximum force required to move the device is 23.5 N

	Avoid positioning the device in an area of direct sunlight exposure.
	Do not position unit near any objects that can generate a strong electrical/magnetic field as this may result in potential electromagnetic interference. The device complies with IEC 60601-1-2.

6.4 Lock/ unlock the castor's brakes

TrueCool has two front and two back castors (see Figure 1). The front castors are fitted with brakes. The brake lever is located over the castors. To lock the castors, firmly press the lever. To release the castors, lift the lever. When the unit is stationary, the brakes must be in the locked position. Release the brakes only while positioning the unit.

	Do not use positioning function of the device while brakes are locked.
	Control unit should be positioned on a flat surface only.

6.5 Storage and transportation conditions

Store TrueCool in a clean and dry area. Disconnect any hoses, blankets, sensors and empty the tank according to subsection 11.2.1. Protect the device to avoid any damage covering it with a protective cover. Secure other accessories. If the device was stored in low temperatures, wait 24 hours to before use. It is recommended to store the device and additional equipment at normal room conditions.

Ambient temperature (Operating)	15-25 °C
Ambient temperature (Storage)	10-40 °C

Ambient temperature (Transportation)	-30-60°C
Humidity (Operating)	30-85%
Humidity (Storage)	30-85%
Humidity (Transportation)	10-85%
Pressure (Operating and Storage)	70-106kPa
Pressure (Transportation)	50-106kPa

6.6 Turning the device on

To connect to the electricity source use the power cord (rear side of the control unit). Switch on TrueCool using the main switch located at the rear of the unit (Figure 4). When the device is switched on, the switch is in the  position and the green light is on. When the device is switched off the switch is in the  position. The main switch is shown below (Figure 15).

	Verify that the voltage switch is set for the local voltage.
	Do not by-pass ground lug. Electrical hazards may result.
	Use an electricity socket with a grounding wire.





Figure 15. Main switch.

Boot – up sequence

When starting TrueCool carries out a boot-up sequence where information about the device is shown. The information that appears on screen must be read by User. If the boot-up sequence was carried out without interruptions, the main screen appears (Figure 22 or 28). If the boot-up sequence was interrupted necessary actions are provided on screen.

 Figure 12. Informative window.	System informs that TrueCool temperature management system must be used in presence of a supervisor.
 Figure 13. Initial tests window.	Initial tests of the condition of the device. IMPORTANT: The initial test sound must be played due to the safety of using.
 Figure 14. Initial test – fault condition (alarm system).	Initial test of alarm system checking has failed.
 Figure 15. Initial test – fault condition (water level).	Initial water level check has failed. The add water informative window appears (figure 19).
 Figure 16. Initial test – fault condition (system check).	Initial alarm system check has failed.

6.7 Filling of the water reservoir

Water should be added when required and before first use. It is recommended to use distilled water. In the event that distilled water is not available then tap water may be used, but this can cause a build-up of lime scale leading to a less effective device in the future.

The water level is checked during the boot up sequence. If the water level is below the limit the device will discontinue boot up and ask to add water.

To add water lift the rear cover (Figure 17), pour water into the inlet whilst watching the onscreen water level marker. Add water to the maximum level. The device will also play an audible sound when the maximum level has been reached.


Figure 17. The water fill opening


Should there be a need to add water during a therapy please do so according to the same above instructions.

The water level marker is shown in Figure 18. When the water level is below the maximum level, the stripes are white instead of blue.

True
Cool™
Temperature management system

Add water until a
beep is heard

Water level

Figure 18. Tank filling – water level marker – boot-up sequence.

	Use distilled water.
	Fill the reservoir only when the control unit is switched on and the water level marker is visible on screen.
	Do not overfill. Overfilling may result in overflow when the water in the blanket drains back into the system when the system is turned off.

6.8.Assembly of hoses and blankets

Hoses

One set of hoses is dedicated for connecting one blanket to the device. Hoses sockets on the device are marked by label with reference number of the hoses (TC1). Figure 8 shows the hose set and its male and female couplings. The Female couplings are used to connect blankets. One male and one female coupling is dedicated to connect hoses to the device (refer to Figure 19a).

Use only hoses recommended by the manufacturer.


Figure 19a. Hoses connection.

Blankets and hoses connection

Hoses are essential to connect blankets to the device. The maximum number of blankets, which may be connected to the device at the same time is 3. The combination of blankets that is used is optional and depends on the actual requirement.

Each blanket has 2 male couplings to be plugged in and 2 female couplings. Figure 19b presents proper blanket and hoses connection. Every blankets is equipped with a valve to prevent water flowing out while the blanket is disconnected. If blankets were properly assembled water will not leak form the connection of blankets and hoses during the work of the device and the therapy.


Figure 19b. Blankets and hoses connection.

Before use read the attached Instruction For Use of the blanket/pad. Follow the manufacturer's instruction in order to safety use blankets with the TrueCool device. Make sure the blanket/pad/mattress is compatible with the device and the hose nozzles fits to the connection socket of the applied part.

7. Patient's preparation

The user must understand all principles of using TrueCool Temperature Management System and all precautions and warnings included in this User Manual must be understood during trainings. It is recommended that user manual as well as instructions for use of each accessory should be reviewed prior to each application.

Before therapy, remove from patient's body (especially from a parts where applied parts to be placed) any liquids, ointments or other agents, as they may influence on the result of the therapy or hazard may occur for the patient. Pay special attention to all transdermal medications as patches or the ointments. Warming may cause too fast drug delivery to the patient's body.

	Warming transdermal medications (patches) can accelerate drug delivery.
	Pay special attention during the course of therapy on patients with ischemic limbs and patients with arterial cross clamping. Long exposure to warming function may lead to burns.
	It is recommended to connect sensors, hoses and blankets when the device is on but before a therapy program has been chosen. To read more about available therapy programs and their selection, refer to chapter 9.

7.1 Sensor application

The available types of the sensors differ in size and method of application. They are two types due to the size: adult and infant, and three types due to the application method: core (to be placed in the patient body), tympanic (to be placed in the ear) and skin (to be placed on the patient body surface); please refer to the section 13.1 for details. As the application method is straightforward, the size should be chosen based on the user experience, and moreover the particular patient, as the patient body may vary and may not be a function of weight, age or height.

Sensors must be connected according to subsection 4.3 to the appropriate sockets. Sensors may be connected while the system is in operation.

Patient temperature sensors may be used in direct connection with patient's body. Used raw materials

for production have been testing according to biocompatibility with human skin and results were positive.

- **Core temperature sensor (esophageal)**

The sensor must be applied accordingly rectal or esophageal. All sensor's head should be applied into the patient. This temperature sensor monitors the internal body temperature and controls the programs based on patient temperature (feedback control). Sensors is of a special type, what ensure the safety of using, high accuracy of measuring temperature and correct work of the sensor. The core temperature sensor is inserted into the rectum and secured with tape to the leg of the patient. Proper applied sensors guarantees safe use for patient and do not cause any irritation of patient's internal tissues. Improper placement can cause a sensor to fall out from the body and it may lead to faulty readings coming from environmental disturbances. Because of the application method there is a risk of flooding the sensor's cable with urine, in that case just dry the cable and continue the therapy.

- **Skin temperature sensor**

The sensor monitors surface patient temperature and must be affixed on patient's skin. The skin temperature sensor usually is taped under the patient's arm on the chest. The usage of skin

temperature sensor is optional and provides only additional information, which is not used in the device's algorithm operation, do not use it in case of any doubts or if the sensor may have a negative impact on the patient condition.

- **Tympanic temperature sensor**

Tympanic temperature sensor provides information about brain temperature and must be inserted in a patient's ear. The tympanic temperature sensor must be applied with carefulness to not irritate or to not damage of patient's eardrum, but its location should ensure appropriate efficiency of work and as accurate measurement of temperature as possible. Make sure to use appropriate size of the sensor to minimize the risk of the potential hazard to the patient. In case of doubts, do not force the sensor to be in place as it may cause damage to the patient eardrum.

The usage of tympanic temperature sensor is optional and provides only additional information, which is not used in the device's algorithm operation, do not use it in case of any doubts, difficulties in properly size selection or if the sensor may have a negative impact on the patient condition.

All type of sensors must be cleaned and disinfected according to the facility's internal procedures before first use. For this purpose, follow advises included in the chapter 11.4 Maintenance of temperature sensors. Single use sensors should be disposed after use according to the manufacturer recommendation in 11.5 Disposal or according to local guidelines. Reusable sensors must be cleaned and disinfected after use according to procedure in 11.4 Maintenance of temperature sensors.

	Do not use any damaged system components, accessories, cables, or sensors. This may lead to inaccurate measurements, hazard to the patient and may damage TrueCool Temperature management System.
	Do not pull the cord to disconnect, always pull from the connector.
	Use only sensors recommended by the manufacturer. Failure to do this will cause incorrect temperature readings and may result in inadequate/inappropriate treatment.
	Single-use sensors are not dedicated to reuse. They have to be disposal after the use. Reuse of the single sensors may lead to faulty measurements, which in the extreme case may cause the death of the patient.
	Check on screen, if sensor(s) are connected properly and the patient's temperature is shown on the screen.
	It is recommended to connect sensors, hoses and blankets when the device is on but no therapy program has been started. To read more about available therapy programs and their selection, refer to chapter 9.
	Working with electronic boards, plugs, and cables requires delicate handling (see Chapter 4.3).

7.1.1. Purpose of sensors

To monitor the patient's temperature, 3 types of sensors are used:

1. Core temperature sensor esophageal for Premium and Lite.

This is the main sensor of the device used to monitor the internal temperature of the patient's body.

2. Tympanic temperature sensor for Premium.

This is an additional sensor that should be placed in the patient's ear and is designed to monitor the temperature of the Tympanic membrane.

3. Skin temperature sensor for Premium.

This is an additional sensor that is used to monitor the patient's skin temperature.

Temperature sensors are available in two types of use: disposable and reusable, and also differ by patient type: adult and for children.

Disposable sensors:

1. Core temperature sensor esophageal for adults	TCP-SCA
2. Core temperature sensor esophageal for children	TCP-SCI
3. Skin temperature sensor for adults	TCP-SSA
4. Skin temperature sensor for children	TCP-SSI

Reusable sensors:

1. Core temperature sensor esophageal for adults	TCP-RCA
2. Core temperature sensor esophageal for children	TCP-RCI
3. Skin temperature sensor for adults	TCP-RSA
4. Skin temperature sensor for children	TCP-RSI
5. Tympanic temperature sensor for adults	TCP-RTA
6. Tympanic temperature sensor for children	TCP-RTI

Caution

- do not pull cord to disconnect (always use the connector);
- do not throw, reusable sensors should be placed gently in the box after use and cleaning;
- do not immerse in liquid;
- storage according to manufacturer's recommendation;
- conduct the cleaning and disinfecting process in keeping with manufacturer instructions;
- handle the sensors gently during cleaning or disinfecting, excessive pressure can stretch the covering and/or break the internal wires.

After the expiration date sensors might not measure with complete accuracy as manufacturer assured. Therefore it is recommended to exchange sensors for new after the expiration date.

Appropriate sensor must be plugged into socket according to colours presented in Figure 10. Tympanic temperature sensor must be plugged in to the yellow tympanic socket (Premium version only), skin temperature sensor to the green socket (Premium version only) and core temperature sensor to the red socket. See Figure 10 to compare the sensors endings and the appropriate sensor sockets. Sensor sockets are located on the left side of the device (Figure 6 – TrueCool Premium, Figure 7 – TrueCool Lite).

rs characteristics:

length for all sensors is: 3,0 m $\pm$ 3,0 cm.

	Factory number (REF)	sensor diameter, mm $\pm$ 1,0 mm	sensor length, mm $\pm$ 1,0 mm	sensor height, mm $\pm$ 1,0 mm	sensor width, mm $\pm$ 1,0 mm
ossible sensors:					
temperature sensor esophageal for adults	TCP-SCA	5,0	18,00	-	-
temperature sensor esophageal for children	TCP-SCI	3,0	12,00	-	-
temperature sensor for adults	TCP-SSA	3,0	20,00	-	30,0
temperature sensor for children	TCP-SSI	3,0	12,00	-	30,0
ossible sensors:					
temperature sensor esophageal for adults	TCP-RCA	5,0	18,0	-	-
temperature sensor esophageal for children	TCP-RCI	3,0	12,0	-	-
temperature sensor for adults	TCP-RSA	9,5	-	4,0	-
temperature sensor for children	TCP-RSI	5,5	-	3,0	-
panic temperature sensor for adults	TCP-RTA	15,0	29,0	-	-
panic temperature sensor for children	TCP-RTI	14,0	22,0	-	-

7.2. Blankets

For manipulations to control the target temperature with the help of the Device for hypothermia with circulating fluid TrueCool, disposable and reusable blankets are used:

TrueCool blankets are designed for temperature therapy in combination with a compatible temperature control device using a water circulation pump with a normal operating pressure in the range of 24 kPa to 35 kPa. Use is limited to temperature regulation under the guidance of a doctor in a medical facility setting.

The disposable blankets, vests, helmets, and TrueCool TCB kits manufactured by EM-MED Sp. z o.o. Sp. K. are designed specifically for TrueCool, and the device achieves full efficiency only with this type of product. It is possible to use these blankets with other hypothermia devices that meet the following criteria:

- a device that generates a pressure equal to or less than 38 kPa,
- adjustable water temperature range not higher than 42° C and not lower than 4° C

It is strongly recommended that you read the user manual that comes with blankets, vests, helmets, and kits.

Disposable blankets, as well as vests, helmets, and TrueCool kits have a connector with a valve or clip to prevent leakage after disconnecting. Blankets, as well as vests, helmets and TrueCool kits supplied by the manufacturer, are intended for use only by one patient and must be disposed of after the procedure in accordance with the internal procedures of the institution.

The blankets, as well as the vests, helmets, and TrueCool kit supplied by EM-MED, can be used in direct contact with the patient's body. The raw materials used for the production were tested for biocompatibility with human skin.

	Before use a blankets or pads read the user manual provided by its manufacturer.
	Do not pull the cord to disconnect, always pull the connector.
	It is recommended to connect sensors, hoses and blankets when the device is on but no therapy program has been started. To read more about available therapy programs and their selection, refer to chapter 9.

Options for disposable blankets, vests, head wraps and sets:

Name	Factory number (REF)	Dimension, cm $\pm 5\%$
Blanket with clamp connections 38 cm x53 cm	TCB-21C	38 cm x53 cm
Blanket with valve connections 38 cm x53 cm	TCB-21V	38 cm x53 cm
Blanket with clamp connections 66 cm x84 cm	TCB-31C	66 cm x84 cm

Blanket with valve connections 66 cm x84 cm	TCB-31V	66 cm x84 cm
Blanket with clamp connections 66 cm x168 cm	TCB-61C	66 cm x168 cm
Blanket with valve connections 66 cm x168 cm	TCB-61V	66 cm x168 cm
TrueCool set includes: vest, head wrap and leg pad	TCS-91V	Head wrap: - width of the top part: 62 cm, width of the lower part: 64,5 cm, height: 45cm; Vest – front width: 51 cm, width at the back: 61 cm, height: 135,5 cm; Leg pad: 66 cm x 84 cm.
Head wrap	TCS-94V	width of the top part: 62 cm, width of the lower part: 64,5 cm, height: 45 cm
Vest	TCS-95V	front width: 51 cm, width at the back: 61 cm, height: 135,5 cm

Reusable TSM blankets are also used:

Reusable Blanket with Valve Connections
1 cm x 145 cm

TCM-51V 51cm x 145cm

Caution

Check the patient's skin condition every 20 minutes (or as directed by a doctor) for individual differences in sensitivity and susceptibility to injury from cold or heat and/or external chemicals or pressure. Patients who are unconscious, who are on long-term therapy, diabetics, children or the disabled, as well as people with sensitive skin areas or poor blood circulation are most at risk.

Precautions and control procedures should be carefully reviewed before each application.

1. Check the blanket for damage.
2. Check the water level in the device, making sure that it is full.
3. Check the bed or operating table for foreign objects and dirt.

4. Unfold the blanket on the bed or operating table or cover the patient.
5. Attach the blanket to the device with hoses.
6. Connect the hoses to the control unit of the device (if it is not already connected).
7. Make sure that there are no folds or kinks in the tube on the blanket.
8. Turn on the device to fill the blanket. The blanket can be filled before or after the patient is on or under it.

Caution

1. Do not use if there are any signs of damage to the blanket.
2. Do not sterilize.
3. Do not operate at temperatures above 42°C or below 4°C.
4. Do not use with a pressure generating device, 38kPa
5. Do not pierce the pads with sharp objects. A puncture will cause water to leak.
6. Do not use pins or other sharp objects to secure the gaskets. Do not allow the gaskets to come in contact with sharp objects.
7. The presence of water or gel on the blanket when the machine is not working can cool the patient.
8. Blankets should be changed if they are dirty.

8. Operating systems screens.

8.1 TrueCool Premium – operating systems screens

8.1.1 Main screen description

After the boot-up sequence has been successfully finished, the main screen appears. The Main screen (Figure 20) is divided into 4 sections: Patient data, Starting therapy, Database and Settings. Patient's data section is used to provide patient's details. This function, as well as using the database, is optional. To pass to next screen (therapy screen) button START must be pressed. On this screen some settings may be provided. Refer to section 9 to read about possible SETTINGS.


Figure 20. Main screen.

8.1.2 Settings screen description

In order to access the settings menu, press the SETTINGS button on the main screen. In the SETTINGS menu (Figure

21) users may adjust the system to their requirements. Possible options are:

- Time/date: select time/date option to set the present time and date.
- Language: select language option to set the language of TrueCool Temperature Management System.

TrueCool

- Maintenance: This option is reserved for authorized service, after selecting this option a PIN access screen appears.
- Producer: select Producer option to read Manufacturer's details.
- Version: select Version option to read Software details.

Figure 21. Settings screen.

TrueCool Premium provides the ability to store patient data that may be subject to local data privacy regulations, and so please follow the appropriate security procedures.

8.1.3 Database description

In order to have an access to the database press the DATABASE button. Database option may be used to download patient's data (sample list shown in Figure 22) or to see a therapy preview by pressing the PLOT button. To export data please PLOT the therapy, insert a flash drive in the USB version and then press export as per on screen instructions. It is recommended to use modern flash drives.


Figure 22. Database screen.

8.1.4 Therapy screen description

The Therapy screen appears after the START button on main screen was pressed. The Therapy screen is divided into graph area, patient data and therapy indication area. The accurate description of individual areas is shown in Figure 23.


Figure 23. Therapy screen description.

The Therapy Screen shows patient data, selected therapy program (Refer to chapter 9 for information on the different therapy programs), measured temperatures (please note that skin and tympanic temperatures will only be displayed if a sensor has been connected), water level, water flow (blue circulating motion).

Measured temperatures are plotted on a graph which can be zoomed in or out using the + and/or - buttons. Left and right arrow buttons are for movements on the time axis when the graph has been zoomed in.

Start button is used to start the therapy and end of therapy button is used to stop the therapy.

Caution: The graphical representations of the patient's physiologic parameters as they are shown by the TrueCool software are tools for displaying the relationships between certain physiologic parameters or the temporal evolution of parameters and do not present any therapy or suggestions for therapy.

8.1.5. Therapy program screen

Therapy Program screen appears after program menu has been selected (Program selection and status area). Three possible program are included into menu: manual water control program, fast automatic program, automatic variable program, temperature change per hour program. The description of therapy programs is enclosed in chapter 9 of this user manual. Figure 24. Program's screen. Figure 24 shows the Program's screen.

Figure 24. Program's screen.

8.1.6. Summary screen

The Summary screen appears after the therapy has been ended. This screen allows the review of the therapy and to analyse the results as well as to export the therapy data. To export the therapy data connect a flash drive press the export button and follow the instructions visible on the screen. It is recommended to use the latest version of flash drive. Press DATABASE button to move on to the Database screen (Figure 22) and plot another therapies results from the list, To return to Therapy Screen press button EXIT. Figure 25 shows the Summary Screen.

Figure 25. Summary screen.

8.2 TrueCool Lite - operating systems screens.

8.2.1 Main screen description

After the boot-up sequence has been successfully finished, the main screen appears. From the Main screen (Figure

26) to start a therapy please select a program from the drop down menu (programs that are available are shown in Figure 27 and are described in chapter 9). Then select a target temperature and press start. When the device has started the water flow indicator will cycle and that status will change to

TrueCool

Page 53 of 68

Rev. 1.4:07.12.2020

Working. Core temperature and water temperature are constantly monitored. To end the therapy press the End Therapy button.


Figure 26. Main screen description – TrueCool Lite.

8.2.2 Therapy program screen

The Therapy Program selection screen appears after program menu has been selected where the user can choose one of three therapy programs: manual water control, fast automatic or automatic variable. The description of each therapy programs is explained in chapter 9.


Figure 27. Therapy Selection Program screen.

9. Description of patient therapy programs

9.1 Manual water temperature control (Premium and Lite)

The program's operation principles are based on controlling the water temperature in the water reservoir. The user defines the target water temperature in a range of 4°C to 42°C using the +/- buttons. The upper and lower limit is fixed and the temperature values cannot be set to the value outside the range. Default value is set at 23°C. The algorithm of the program, based on the readings from the water temperature sensor located in the water reservoir will decide what function (warming or cooling) must be used to achieve the desired water temperature (TSET). Depending

TrueCool

on the TSET in relation to water temperature, the appropriate operation mode will be indicated on the device's screen. When the water temperature is in the defined range of $TSET \pm 0.5^{\circ}\text{C}$, the algorithm will tend to keep the water temperature within stated tolerance by turning off and switching on appropriate function. The algorithm is constantly monitor the water temperature until therapy is stopped.

	Monitor patient temperature constantly.
	Do not leave TrueCool Temperature Management System unattended while the therapy is being carried out.
	Make sure that the vital functions of the patient are always monitored by independent medical equipment while the therapy is carried out.

9.2 Fast automatic (Premium and Lite)

The program's operation principles are based on feedback control from the patient core temperature sensor. The user defines the target patient temperature. The target patient temperature is set with the +/- buttons. The Target patient temperature may be defined by the user in the range of: $30^{\circ}\text{C} - 40^{\circ}\text{C}$. The upper and lower limit is fixed and the temperature values cannot be set to the value outside the range. Default value is set at 36°C . The algorithm of the program, based on the readings from the patient core temperature sensor will decide what function (warming or cooling) must be used to achieve desired patient core temperature (TSET). Depending on the patient core temperature in relation to TSET, the appropriate operation mode will be indicated on the device's screen. When the patient core temperature is in the defined range ($TSET \pm 0.05^{\circ}\text{C}$ for TrueCool Premium, and $TSET \pm 0.1^{\circ}\text{C}$ for TrueCool Lite), the algorithm will tend to keep the patient core temperature within stated tolerance by turning off and switching on appropriate function. In this program, the water temperature which can be achieved by the system is limited to the range of 4°C to 42°C (minimum and maximum value to be achieved by the device in normal use). If the actual function of the device is warming, then the water tends to reach the maximum available value of temperature (42°C), if the actual function of the device is cooling, then the water tends to reach the minimum available value of temperature (4°C). If the water reaches the minimum or maximum value of temperature, the algorithm will tend to keep it at constant level, within stated tolerance (as it is described in Chapter 9.1). The algorithm is constantly monitor the patient and water temperature until therapy is stopped.

Figure 28 shows the program selection window for TrueCool Premium and Figure 29, program selection window for TrueCool Lite.

	Make sure that the core temperature sensor is connected before commencing therapy.
	Always check whether the core temperature is visible in temperature indication area.
	Monitor patient temperature constantly.


Figure 28. Fast automatic program – TrueCool Premium.



Figure 29. Fast automatic program – TrueCool Lite.

9.3 Automatic variable gradient (Premium and Lite)

The program's operation principles are based on feedback control from the patient core temperature sensor. The user defines the target patient temperature. The target patient temperature is set with the +/- buttons. The Target Patient temperature may be defined by the user in the range of: 30°C – 40°C. The upper and lower limit is fixed and the temperature values cannot be set to values outside the range. Default value is set at 36°C. The algorithm of the program, based on the readings from the patient core temperature sensor will decide what function (warming or cooling) must be used to achieve desired patient core temperature (TSET). Depending on the patient core temperature in relation to TSET, the appropriate operation mode will be indicated on the device's screen. In this program, the water temperature which can be achieved by the system is limited to the "boundary" range of 4°C to 42°C (minimum and maximum value to be achieved by the device in normal use) and additionally by the user defined gradient of the temperature in the range of 0°C to 15°C (default value of the gradient is set at 10°C).

In other words, gradient temperature additionally limits the water temperature. If the actual function of the device is warming, then the water tends to reach the maximum available value of temperature which is $TPATIENT + TGRADIENT$. The sum of the $TPATIENT + TGRADIENT$ cannot be higher than 42°C, if so, the water temperature is kept at the level of 42°C. If the actual function of the device is cooling, then the water tends to reach the minimum available value of temperature which is $TPATIENT - TGRADIENT$. The difference of the $TPATIENT - TGRADIENT$ cannot be lower than 4°C, if so, the water temperature is kept at the level of 4°C. If the water reaches the minimum or maximum allowable value of temperature, the algorithm will tend to keep it at constant level, within stated tolerance (as it is described in Chapter 9.1). The above conditions applies to a stage, when the desired patient temperature is not reached since the therapy is started. When the patient core temperature is in the defined range ($TSET \pm 0.05^\circ\text{C}$ for TrueCool Premium, and $TSET \pm 0.1^\circ\text{C}$ for TrueCool Lite), the algorithm will tend to keep the patient core temperature within stated tolerance by turning off and switching on appropriate function, but since the patient temperature reached exactly TSET the water temperature is not limited anymore by the gradient value but only by "boundary" conditions as in Fast automatic program (4°C to 42°C). The new conditions remain unchanged till the program is stopped or its parameters are changed. The algorithm constantly monitor the patient and water temperature until therapy is stopped.

Figure 30 shows the program selection window for TrueCool Premium and Figure 31, program selection window for TrueCool Lite.

	Make sure that the core temperature sensor is connected before commencing therapy.
	Always check whether the core temperature is visible in temperature indication area.
	Monitor patient temperature constantly.

9.3.1 Automatic variable gradient + delayed boost

The automatic variable gradient therapy program has a delayed boost option available for the user (marked by default). This option allows the algorithm to increase the primary value of gradient by 5°C, if the patient core temperature has not been changed by 0.1°C within 30 minutes since the option has been marked.

Since the option is marked, the values of the patient temperature are recorded with one minute interval, which means that the first check of a condition is performed 29 minutes (30 available values of patient core temperatures) from the time when the option is marked. If the patient core temperature has been changed by 0.1°C, the primary value of the gradient remains unchanged, and condition is check every sixty seconds till the option is marked or therapy stopped. If the patient core temperature has not been changed by 0.1°C, the primary value of the gradient is then modified by increase its value by 5°C, since the value has been changed, all the previously saved values of patient core temperature are deleted and the cycle repeats like the boost option would be marked for the first time (we need to record a new 30 values of the patient core temperature).

Figure 32 shows the program selection window for TrueCool Premium and Figure 33, program selection window for TrueCool Lite.


Figure 32. Automatic variable gradient program delayed boost option activated – TrueCool Premium.


Figure 33. Automatic variable gradient program delayed boost option activated – TrueCool Lite.

	Make sure that the core temperature sensor is connected before commencing therapy.
	Always check whether the core temperature is visible in temperature indication area.
	Monitor patient temperature constantly.

9.4 Temperature change per hour (Premium only)

The program's operation principles are based on feedback control from the patient core temperature sensor. The user defines the target patient temperature. The target patient temperature is set with the +/- buttons. The Target patient temperature may be defined by the user in the range of: 30°C – 40°C. The upper and lower limit is fixed and the temperature values cannot be set to the value outside the range. Default value is set at 36°C. The Target patient temperature in this program must be reached in a specific manner – by increasing or decreasing the patient core temperature every 60 minutes, by user defined value from the range of 0.05°C – 0.5°C. The general idea of this program can be described, as the fast automatic program, in which every 60 minutes the system automatically, increase or decrease the desired patient temperature (T^{HSET}) by user-defined value, starting from the actual reading of the patient core temperature. Within each hour, the patient core temperature is kept at the constant level within specified range using appropriate function (warming or cooling). Depending on the patient core temperature in relation to TSET, the appropriate operation mode will be indicated on the device's screen. The changing of the "local" desired patient core temperature (T^{HSET}) will be performed until the "final" desired patient core temperature (TSET) is reached ($T^{HSET}=TSET$). In this program, the water temperature which can be achieved by the system is limited to the range of 4°C to 42°C (minimum and maximum value to be achieved by the device in normal use). If the actual function of the device is warming, then the water tends to reach the maximum available value of temperature (42°C), if the actual function of the device is cooling, then the water tends to reach the minimum available value of temperature (4°C). If the water reaches the minimum or maximum value of temperature, the algorithm will tend to keep it at constant level, within stated tolerance (as it is described in Chapter 9.1). The algorithm is constantly monitoring the patient and water temperature until therapy is stopped.

Figure 34 shows the program selection window.


Figure 34. Temperature change per hour program.

▲	Make sure that the core temperature sensor is connected before commencing therapy.
▲	Always check whether the core temperature is visible in temperature indication area.

9.5 Therapy ending

9.5.1 TrueCool LITE

When the End Therapy button is pressed the below screen will appear. Refer to chapter 6.


Figure 35. End of therapy screen.

9.2 TrueCool PREMIUM

When the End Therapy button is pressed the main menu screen will appear.






Figure 36. Main menu screen.


10. Alarms and information windows

Alarms are described in section 5.5.

Information windows

Window	Description
End of therapy confirmation window 	End of therapy confirmation window contain two action buttons – exit/change button "<" and confirmation button "OK".
No USB window 	If the Flash drive icon is pressed with no USB drive connected to the USB Socket, then the error will be displayed. It is necessary to press the action button "Back" to close the window.
Save therapy window 	If the Flash drive icon is pressed (see point 32) with USB drive connected to the USB Socket, then the confirmation window is displayed. To accept data export/saving to the external drive, the "Save" button needs to be pressed. It is necessary to press the action button "Back" to close the window without saving the data to the USB drive.
Saving succesful window 	If the export/saving data is complete, the confirmation window is displayed. It is necessary to press the action button "Back" to close the window.
No core Sensor Window	The windows is displayed if no patient core temperature sensor is connected to the device and the operator wants to start a therapy program dependent on patient temperature.

Window	Description
--------	-------------

Symbol  is used to indicate that action is required.

Information icons

Window	Description
Maintenance required 	Indication will show up next to the SETTINGS button, if the moto-hours counter reach the appropriate value.
UV replacement required 	Indication will show up next to the SETTINGS button, if the UV hours counter reach the appropriate value.

11. General user maintenance/cleaning principles

This section describes the general requirements for maintenance of the unit which can be performed by the user to keep the system in a good working condition between annual inspections (Technical Service Inspection, which should be carried out minimum once a year). The TrueCool housing is constructed of painted, resin polyester-glass laminate that should be cleaned quarterly with water and a mild detergent solution.

The device will also show the service reminder icon indicating that the general maintenance is required due to the device or UV LEDs running time. The appropriate symbol will be displayed next to the SETTINGS button on the main screen.

	Always unplug the unit before accessing internal components during service. Failure to unplug the unit could result in electric shock.
	The repair, calibration, and servicing of the TrueCool should be performed by qualified Medical Equipment Service Technicians, Certified Biomedical Electronics Technicians, or Certified Clinical Engineers familiar with good repair practices for servicing medical devices, and in accordance with instructions contained in this manual and the Service Manual. Improper repair can result in damage to the TrueCool system and patient injury.

11.1 Cleaning and disinfection

Cleaning and disinfection of the device.

Clean the exterior body of the unit, fluid delivery hoses and power cord using a soft cloth and mild detergent or disinfectant according to hospital protocol. Always use conventional hospital-approved disinfectants that do not contain alcohol. Thoroughly wipe off any residue of the cleaning solutions. Special attention must be paid to the ventilation holes which should be free from dust or any other dirt which may block the airflow.

Cleaning and disinfection of blankets.

Disposable blankets are not subject to cleaning or disinfection. If they get dirty in any way, it is recommended to throw away the blanket and replace it with a new one.

Before the first use, the reusable blanket can be cleaned with a damp cloth with water and, for example, with soap. Disinfection of a contaminated area on a reusable blanket can be carried out using a variety of available means that are effective against a wide range of bacteria and microorganisms. Products based on sodium hypochloride are safe to use. If you get dirty, collect as many pollutants as possible and throw them away. Using germicidal wipes or a bactericidal spray with disposable towels wipe the affected area. Depending on the number and area of the lesion, it may take several attempts to thoroughly clean the blanket. After the area has been wiped with a disinfectant, allow the air to dry.

Recommended solutions / detergents:

- Products based on sodium hypochloride
- disposable disinfectant / sanitizer wipes CaviWipes1 line ®
- any product of the CaviCidel ® line
- any product of the Sani-Cloth ® line

Cleaning and disinfection of sensors.

The sensors must be clean before use.

Follow the procedure:

1. Use a damp cloth with water and neutral detergent. Clean all exposed surfaces of the cable and sensor.
2. Use another wet cloth with sterile or distilled water to remove all chemicals.
3. After cleaning, it is recommended to dry the sensor and cable.

WARNING!

Never immerse the temperature sensor in a liquid

DISINFECTION

The preferred method is ethylene oxide. You can also use 3% hydrogen peroxide, 70% isopropanol, 70% ethanol, CIDEX, a 1:10 solution of sodium hypochlorite.

- any product of the Sani-Cloth ® line

11.2 Maintenance of the water reservoir

11.2.1 Draining the reservoir

It is recommended to replace water in the reservoir every quarter with distilled water. To drain the reservoir place the device under a sanitary drain or place a container under the valve (Figure 1). Open the valve and wait until container is full. Close the valve, empty the container and repeat the procedure until no more water flows out of the reservoir. Proceed to add water according to chapter 6.7.

	Always drain water or dispose of to a sanitary drain because bio-contaminants may be present in the unit's water supply
	Remember to clean/disinfect the control unit/hoses/sensors after each finished therapy.

11.2.2 Fluid circuit disinfection

11.2.2.1 Quarterly

It is recommended to conduct a quarterly disinfection by switching the unit on and letting it run for 15 minutes in manual mode at water set temperature 22°C.

It is not recommended to use de-ionized water or alcohol as it may damage the internal parts of the device.

The UV disinfection will clean the water in the reservoir.

UV disinfection will purify the water in the tank.

The vents at the back of the unit should be vacuumed during quarterly maintenance to keep the unit well ventilated during operation.

11.2.2.2 Yearly

The procedure requires the use of one of the substances listed in this manual (see Annex A at the end of manual) or alternatively bleach may be used.

It is recommended, to use a single use blanket/pad/mattress during the procedure, so the hoses can be also flushed with water. This will also improve the flushing efficiency through the water circulation system.

The FLUID CIRCUIT DISINFECTION procedure is NOT to be performed during patient therapy.

Follow the instructions:

1. Turn the device off.
2. Drain water out of the system.
3. Fill the reservoir with the appropriate amount of water and chemical compound from the manual. It is recommended to use around 300 ml of drug (or alternatively 300 ml of bleach) for full water reservoir (around 10 liters).
4. Turn on the device and set the Manual water temperature control program to 35 – 40°Celsius.
5. Press the START button and leave the device for a few minutes in operation with the compound in the system.
6. If the blanket/pad/mattress is to be used during procedure, attached it to the side couplings and leave it for another several minutes.
7. Drain the reservoir.
8. Repeat the steps 3-7 three times (before repeating the steps, change the hose set if more than one is available).
9. After the third rinse, drain the unit and fill the tank with distilled or tap water and set the program.
10. Drain the water out of the system.
11. Check the poured water with pH strips or other appropriate method for used chemical compound detection.
12. If the tests detect chemical compound, which was used, repeat the steps 9-11.
13. Refill unit with distilled water.

	For safe handling and use of chemicals follow manufacturer guidelines.
	Users should not use cleaning or decontamination methods different from those recommended by EM-MED without first checking with EM-MED that the proposed methods will not damage the equipment.
	Do not use strong chlorine solutions.

11.3 Maintenance of hoses

The hoses do not require a special maintenance, however, it is recommended to flush it with water during yearly maintenance as described in subsection 11.2.2.2.

11.4 Maintenance of temperature sensors

All types of sensors should be clean before first use with a damp cloth and mild detergent. When not in use, the sensor(s) should be loosely coiled or placed back in the original box and stored at room temperature.

After use reusable sensors should be disinfected, for this purpose use isopropyl alcohol. Damp cloth with recommended agents and disinfect gently the applied part of the sensor, wire and then plug connector (excluded metal part of the connector). Pay attention to thoroughly disinfect sensor's cavities, especially in case of tympanic sensors with an irregular shape. Not proper or not effective disinfection process of reusable sensors may be cause of infection of the patient, thus is strongly recommended to follow this instruction and use proven agents.

Avoid sensors contact with strong aromatic chlorinated, ketone, ether or ester solvents, additional further avoid soaking the wire in alcohol, harsh detergents or highly alkaline substances. The application of ether or acetone causes the loss of flexibility in the cable.

Remember, while cleaning or disinfecting, the sensor should be handled gently. Excessive pressure can stretch the covering and/or break the internal wires.

	Do not autoclave.
	Do not immerse the temperature sensor jack connector in liquid.
	Do not immerse temperature sensors or connecting cables in liquid.
	Working with electronic boards, plugs, and cables requires delicate handling (see Chapter 4.3).

11.5 Disposal

Accessories of the medical device that have come in contact with patients contain the risk of bio-contamination. In purpose to disposal of patient temperature sensors, single-use and reusable after of their useful lives, follow the local State and hospital guidelines regarding to disposal of the medical device.

For disposal blankets please follow instructions in the their user manual.

The TrueCool has cooling agent R290, which may lead to biological contamination of the environment and its incorrect disposal may lead to fire explosion, therefore before the disposal process of the medical device the R290 must be removed only by trained service technician. Additional the TrueCool consists of electronic components, which also may be hazardous for the environment. The device should be disposed of in accordance with the local State and hospital guidelines regarding to disposal of the medical device. In other case the medical device may either be returned to the supplier for recycling.

Remove from the device the cooling agent R290 before disposal. Action must be conducted only by trained service technician. Incorrect disposal may lead to fire or explosion.

12. Technical specification

Table 10. Physical characteristic

Physical characteristic	TrueCool Premium	TrueCool Lite
	10" display with touch panel	7" display with touch panel
Construction	Painted, polyester-glass resin laminate housing. Safety glass on the front panel. Divided water reservoir (approximately 10 liters). Bottom air discharge. Big handle at the front of the device. Base of the housing as a bumper guards. Hook for the hoses at the rear wall of the device. Four 75 mm conductive swivel-castors, two (front) with blocking mechanism. Stainless steel frame construction.	
Connections	3 temperature probes sockets (YSI 400 Series Probes) 1 USB socket 3 inlet and 3 outlet couplings for hoses connection (quick disconnect couplings).	1 temperature probe socket (YSI 400 Series Probes) 1 USB socket 3 inlet and 3 outlet couplings for hoses connection (quick disconnect couplings).
External Dimensions	Width: 50 cm Depth: 55 cm; with hook Open: 70 cm Height: 102 cm	
Weight (device)	66 kg	
Maximum force required to move the device	23,5N	

Table 11. Use conditions.

Use conditions	TrueCool Premium	TrueCool Lite
Ambient temperature (Operating)	15-25 °C	
Ambient temperature (Storage)	10-40 °C	
Ambient temperature (Transportation)	-30-60°C	
Humidity (Operating)	30-85%	
Humidity (Storage)	30-85%	
Humidity (Transportation)	10-85%	
Pressure (Operating and Storage)	70-106kPa	
Pressure (Transportation)	50-106kPa	

Table 12. Characteristic.

Characteristic	TrueCool Premium	TrueCool Lite
Type of protection against electric shock	Class 1	Class 1
Operating mode	lasting mode	lasting mode
Electricity Input Power	220 – 240V, 50Hz	220 – 240V, 50Hz
Max Power Input Total	1.5kVA	1.5kVA
Input rating	1470VA	1470VA
Maximum Current	6.5A	6.5A
On-board FUSES	16A, 250VAC, time lag; 632.330	16A, 250VAC, time lag; 632.330
Degree of protection against dangerous penetration of water or solid particles	IP2X	IP2X
Compressor NEK6213U	P=471W; In=3,47A	P=471W; In=3,47A
Condenser fan	P=38W	P=38W
Valve coil	P=10W	P=10W
Heater	P=800W	P=800W
Pump	P=60W; In=0,24A	P=60W; In=0,24A
Blanket normal operating pressure	Up to 6 PSI	
Disinfection system – biological safety	Build-in UV LEDs diodes in a closed circulatory system ensuring no photo biological hazard.	
Refrigerant (chladogen)	R290 (C3F8)	
Coolant mass	135g	
Water temperature range (available to be set by the operator)	4° to 42°	
Patient temperature range (available to be set by the operator)	30° to 40°	
Primary Overheating Alarm	For water temperature > 43.0°C±0.5°C	
Secondary Overheating Alarm	For water temperature ≥44.0°C±0.5°C	
Primary Overcooling Alarm	For water temperature < 3.0°C±0.5°C	
Secondary Overcooling Alarm	For water temperature ≤ 2.0°C±0.5°C	
Accuracy of the Patient core temperature sensor	0.05°C	0.1°C
Accuracy of the water temperature sensor	0.5°C	0.5°C
Resolution of temperature values	0.01°C	0.1°C
Device Classification	IIB	
Degrees of protection against electrical shock	BF	

Table 13. External dimensions

External dimensions	Characteristic
Device TrueCool Premium and Lite	Table 10. Physical characteristic.
Power cable length	2 m
Sensors - cable length	3m
Hoses length	3m

Table 14. Optional functions.

Optional CFM/EEG device requirements – Premium version only	
Regulatory requirements	The device must have CE sign.
	The device must be a medical device with EC certificate.
Technical requirements	The device must be fitted with USB socket (receiving signal from TrueCool device).
	The weight of the screen mounted on the support arms cannot exceed 15 kg.
	The screen must have a proper fixation holder.

13. List of accessories

List of accessories depends on the TrueCool model and customer order.

- Temperatures sensors: core temperature sensor, tympanic temperature sensor, skin temperature sensor, single use or reusable (excluding tympanic).
- Hoses up to 3 sets.
- Holder for the CFM/EEG device.

14. Electromagnetic interference

For safe use of the TrueCool it is required to keep TrueCool at safe distance from devices emitting radio frequency. Interruptions in power supply affect the functionality of the system and can damage the device.

Make sure to read the messages in order to ensure correct reactivation of the machine.

15. LED UV emitters

LED UV Emitter is a revolutionary, energy efficient light source for sterilization and medical application, combining the lifetime and reliability advantages of Light Emitting Diodes with the brightness of conventional light source.

A UV disinfection system is installed in the dual water reservoir which works in an automatic 4-hour a day cycle irrespective of a device to patient connection, hindering pathogen growth. The UV disinfection system has been tested and its effectiveness has been confirmed by tests performed by an independent laboratory.

UV LEDs operation time

UV operation time counter (UV counter) is turned on only if UV diodes are turned ON. UV counter is used to check the time, when the UV diodes boards should be replaced with new ones. If the UV counter (Counter (UV LEDs)) is equal or higher than 700 hours, the appropriate indication is displayed on the screen

Indication will show up next to the SETTINGS button, if the UV hours counter reach the appropriate value.

The RESET button is used to reset the value of the Counter (UV LEDs) to zero any time the reset button is pressed by a service technician.

We recommend to keep the minimum spares stock.

Before taking out the UV LED modules from the quartz tubes mounted into the water reservoirs, make

sure the modules are disconnect from the module PCB as it may lead to UVC light exposure.
Direct exposure of skin and eyes to UVC radiation from the UVC diode may cause painful eye injury and burn-like skin reactions. Never look directly at a UVC light source, even briefly.
In Each Premium and Lite model uses 2 UV LED modules.

*Features

- I. C. compatible
- RoHS compliance and Pb free
- Lower operating costs
- Reduced maintenance costs

Available part numbers

Nominal peak wavelength (nm)	Forward Voltage
Wp 280	6.2 V

Absolute Maximum Rating at Ta=25°C

Parameter	Symbol	Rating	Unit
Power Dissipation	Po	0.3	W
DC Forward Current	IF	40mA	mA
Operating Temperature Range	Topr	-30~+80	°C
Storage Temperature Range	Tstg	-40~+100	°C

Notes: Operating the LED under reverse bias condition long time might result in damage or failure of the components

Electro-Optical Characteristics at Ta=25°C

Parameter	Symbol	Values		Test Condition	Unit
Forward Voltage	Vf	Max.	7.5	If = 20mA	V
		Typ.	6.2		
		Min.	4.4		
Radiant Flux	Φe	Max.	3.8	If = 20mA	mW
		Typ.	2.4		
		Min.	1.4		
Peak Wavelength	Wp	Max.	290	If = 20mA	nm
		Min.	270		
Viewing Angle	2Θ 1/2	Typ.	125	If = 20mA	°
Thermal Resistance (Junction-scider point)	Rth j-s	Typ.	35	If = 20mA	K/W

Notes:

- Radiant flux (Φe) is the total Radiant flux output as measured with an integrating sphere.
- Forward voltage (Vf) measurement tolerance is ± 0.1 V
- Peak Wavelength measurement tolerance is ± 3 nm
- Forward Radiant Flux measurement tolerance is ± 10%

Bin Code list

Vf: Forward Voltage Spec. Table		
V _F Bin	V _f (V) @ I _f = 20mA	
	Min.	Max.
VY	4.4	5.1
V0	5.1	5.8
V1	5.8	6.5
V2	6.5	7.0
V3	7.0	7.5

Tolerance on each Forward Voltage bin is +/-0.1V

Φ _e : Radiant Flux Spec. Table		
Φ _e Bin	Φ _e (mW) @ I _f = 20mA	
	Min.	Max.
X1	1.4	-

Tolerance on each Radiant Flux bin is +/-7%

W _p : Peak Wavelength Spec table		
W _p bin	W _p (nm) @ I _f = 20mA	
	Min.	Max.
W1	270	290

Tolerance on each Peak Wavelength bin is +/-3nm

Notes:

1. Bin classification code is marked on each packing bag.

Recommend Printed Circuit Board Attachment Pad**Notes:**

1. The LEDs can be soldered using the reflow soldering or hand soldering method. The recommended hand soldering condition is 300 °C max. and 2 secs max/ only once and the operation of reflow soldering is only up to three times at maximum.

Notes:

1. Empty component packets sealed with top cover tape.
2. 7 inch reel-maximum 1000 pieces per reel. (Min. packing quantities are 100 pieces for remainders)
3. The maximum number of consecutive missing lamps is two.
4. In accordance with EIA-481-1-B specifications.


CAUTIONS

Cleaning


Use alcohol-based cleaning solvents such as isopropyl alcohol to clean the LED if necessary. Do not use unspecified chemical liquid to clean LED, it could harm the package.

Drive Method

An LED is a current-operated device. In order to ensure intensity uniformity on multiple LEDs connected in parallel in an application, it is recommended that a current limiting resistor be incorporated in the drive circuit, in series with each LED as show in Circuit below

Circuit model A

Circuit model B

(A) Recommended circuit

(B) The brightness of each LED might appear different due to the differences in the I-V characteristics of those LEDs.

© This LED should be used under forward current, particularly not used under continuous reverse current to avoid the damage to LED

The electrode pads are plated with gold, but it is still not recommended to the use under any of the following conditions, please confirm the performance and reliability are well enough if you use it under any of the following conditions.

- Do not use sulfur-containing materials in commercial products including the materials such as seals and adhesives that may contain sulfur.
- Do not put this product in a place with a lot of moisture (over 85% relative humidity), dew condensation, briny air, and corrosive gas (CL, H₂S, NH₃, SO₂, NOX, etc.).

Manual Handling Remark

The LED should only be picked up by making contact with the sides of LED body. It should not put any pressure on the lens either by finger or any hand tool. Do not puncture or push the lens.

ESD (Electrostatic Discharge)

- Static Electricity or power surge will damage the LED. Suggestions to prevent ESD damage:
 - Use a conductive wrist band or anti electrostatic glove when handling these LEDs.
 - All devices, equipment, and machinery must be properly grounded.
 - Work tables, storage racks, etc. should be properly grounded.
 - Use ion blower to neutralize the static charge which might have built up on surface of the

LED's plastic lens as a result of friction between LEDs during storage and handling.

ESD-damaged LEDs will exhibit abnormal characteristics such as high reverse leakage current, forward voltage, or "no light up" at low currents. To verify for ESD damage, check for "light up" and V_F of suspect LEDs at low currents.

Training and Certification

1. Working area is ESD-certified.
2. Training records and re certification dates monitored.

Static-Safe Workstation & Work Areas

1. Static-safe working stations or work-areas have ESD signs.
2. All surfaces and objects at all static-safe workstation and within 1 ft measure less than 100V.
3. All ionizer activated, positioned towards the units.
4. Each work surface mats grounding is good.

Personnel Grounding

1. Every person (including visitors) handling ESD sensitive (ESDS) items wear wrist strap, heel strap, conductive shoes with conductive flooring.
2. If conductive footwear used, conductive flooring also present.
3. Garments, hairs or anything closer than 1 ft to ESD items measure less than 100V.
4. The wrist strap or heel strap/conductive shoes are checked daily and result recorded.
5. All wrist strap or heel strap checkers calibration up to date.

Device Handling

1. Each ESDS items identified by EIA-471 labels on item or packaging.
2. No static charge generators (e.g. plastics) inside shielding containers with ESDS items.
3. All flexible conductive and dissipative package materials are inspected before reuse or recycles

Storage

This product is qualified as moisture sensitive Level 3 per JEDEC J-STD-020 Precaution when handling this

moisture sensitive product is important to ensure the reliability of the product.

The package is sealed:

The LEDs should be stored at 30°C or less and 85%RH or less. And the LEDs are limited to use within one year, while the LEDs is packed in moisture-proof package with the desiccants inside.

The package is opened:

The LEDs should be stored at 30°C or less and 60%RH or less. The LEDs are limited to solder process within

168hrs. we recommended to baking LEDs at 60°C at least 24hours before solder. To seal the remainder LEDs

return to package, it's recommended to be with workable desiccants in original package.

Safety

DO NOT directly look at the UV light from LED or optical instrument because it is harmful to human eyes. And it is recommended to use UV goggles protective glasses for avoiding eye damage.

1.6. Shelf life

Device: The service life of the TrueCool circulating fluid hypothermia apparatus is 5 years from the date of manufacture, with proper care and attention.

Temperature sensors: The service life of each of the patient's temperature sensors is 3 years from the date of manufacture, with proper care and attention. After the expiration date, the sensors may not measure accurately. It is recommended to replace the sensors with new ones after the expiration date.

Blankets, TrueCool set includes: vest, head wrap and leg pad, head wrap, vests: The service life of each of the blankets, sets of True Cool in the composition: helmet for the head, vest, blanket for the feet with valve connections, helmets, vests is 5 years from the date of manufacture with proper care and attention. After the expiration date, the sensors may not measure accurately. It is recommended to replace the sensors with new ones after the expiration date.


The warranty period for the device, sensors and blanket, as well as the TrueCool set, is 12 months from the date of receipt of the product by the customer.

Symbols used on the labeling

The symbols used in the labels of the TrueCool circulating fluid hypothermia device with accessories, Premium and Lite models, as well as sensors, blankets, and True Cool sets are listed below:

Manufacturer

Date of manufacture

Catalog number

Batch number

Not for general disposal (separate collection of electrical and electronic equipment).

See user manual and instructions for using.

Warning (CAUTION)

Type BF Equipment

Refer to user manual/read the instruction

The product complies with the requirements of Directive 93/42 / EEC

IP2X Device is rated to a degree of Ingress Protection (IP) 2X

Keep dry

Fragile, handle with care

Keep away from sunlight

Single use


Do not sterilize

20 min. Check the patient's skin condition every 20 minutes

Use before

Temperature limit

SN

Serial number

Does not contain natural latex

Medical Devices

Refrigerant

Chladogen R290

Fire hazard. Highly flammable substances

Do not dispose of together with household waste

Packaging

Device: The device is packed in a non-sterile form in a separate cardboard box made of 5-layer cardboard with the appropriate label, which contains the product name, batch number. On all sides, the device is protected by protective tabs made of polyethylene. Packing Units: 1 unit in a box with cable and user manual.

Sensors:

Disposable sensors are packed in a plastic package with a Zip lock. Packing units: 1 unit per pack

Reusable sensors are packed in cardboard packaging. Packing units: 1 unit per pack

Blankets:

The disposable blanket is packed in a plastic bag. Packing units: 1 unit per pack.

The package contains an instruction manual with the name, article number, batch number and expiration date. Next, the blankets in the amount of 10 units are packed in a cardboard box with a label containing also the name, article, batch number and expiration date.

Blankets, TrueCool set includes: vest, head wrap and leg pad, head wrap, vests:

Each set is packed in a cardboard box with enclosed operating instructions with the name, article number, batch number and expiration date. Packing units: 1 set per pack. Further, the sets in the amount of 4 units are packed in a cardboard package with a label containing also the name, article, batch number and expiration date.

Head wrap:

The head wrap for the head is packed in plastic packaging. Packing units: 1 unit per pack. The package contains an instruction manual with the name, article number, batch number and expiration date. Further, the helmets in the amount of 10 units are packed in a cardboard box with a label containing also the name, article, batch number and expiration date.

Vest:

The vest is packed in a plastic bag. Packing units: 1 unit per pack. The package contains an instruction manual with the name, article number, batch number and expiration date. Further, the vests in the amount of 5 units are packed in a cardboard box with a label containing also the name, article, batch number and expiration date.

Reusable Blanket:

The reusable blanket is packed in a plastic bag with enclosed operating instructions with the name, article number, batch number and expiration date.

17. List of applicable standards

Standart	Name
EN 60601-1:2006 / A1:2013	Medical electrical equipment - Part 1: General requirements for basic safety and essential performance
EN 60601-1-2:2015	Medical electrical equipment. General requirements for basic safety and essential performance. Collateral Standard. Electromagnetic disturbances. Requirements and tests


EN 60601-1-6:2010	Medical electrical equipment - Part 1-6: General requirements for basic safety and essential performance - Collateral standard: Usability
EN 60601-1-8: 2007 / A11: 2017	Medical electrical equipment - Part 1-8: General requirements for basic safety and essential performance - Collateral Standard: General requirements, tests and guidance for alarm systems in medical electrical equipment and medical electrical systems
EN 60601-1-10:2008	Medical electrical equipment - Part 1-10: General requirements for basic safety and essential performance - Collateral standard: Requirements for the development of physiologic closed-loop controllers
EN 62304: 2006	Medical device software - Software life-cycle processes
EN 62366: 2008	Medical devices - Part 1: Application of usability engineering to medical devices
EN 1041:2008	Information supplied by the manufacturer of medical devices
EN ISO 14971:2012	Medical devices - Application of risk management to medical devices
EN ISO 10993-1:2009	Biological evaluation of medical devices—Part 1: Evaluation and testing within a risk management process
EN ISO 15223-1: 2016	Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements

Other applied standards

Standard	Rationale
EN 14971:2019	Medical devices – Application of risk management to Medical devices
EN IEC 63000:2018	Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances
EN 62353:2014	Medical electrical equipment – Recurrent test and test after repair of medical electrical equipment
EN 14155:2011	Clinical investigation of medical devices for human subjects — Good clinical practice

Exclusions

Standard	Rationale
EN 62471:2008 Photobiological safety of lamps and lamp systems	The standard concerns photo biological safety of use lamps and lamp systems together with lighting fittings and how they effects on people (retinal damage, which can occur in a photochemical or thermal way). TrueCool is not a device used to illuminate anything. The LED operation is only auxiliary and is not the essence of the medical device, therefore the norm is not apply. The UV LEDs have been used as an additional safety system in the event of water contamination in a closed circulation system. Biological safety is ensured by the tightness of the water circulation system. No single

	person is exposed to the radiation. The LEDs are in a closed system so there is no photo biological hazard.
Machinery Directive 2006/42/EC	According to definition, a machine is to be understood as a device equipped with a propulsion mechanism, consisting of parts which are coupled, of which at least one part carries out the motion, joined together for a specific application. TrueCool is stationary product, and it doesn't fulfill definition of machine acc. to 2006/42/EC. There are no moving parts inside it. Outside wheels are used only for positioning the device at the treatment room.
EN 60601-2-26:2003	The standard concerns particular requirements for the basic safety and essential performance of electroencephalographs. The manufacturer is not supplier of CMF/EEG. The TrueCool on custom made could be equipped only with the CMF/EEG holder. Moreover manufacturer informs in User Manual if customer would like to add CMF/EEG, the assembled device must be provide with CE mark.
EN 80601-2-35:2009	The primary function of the device TrueCool device is the cooling function which is assisted by the heating warming function to reduce the intensity / speed of the patient's cooling. It is true that the device enables the regulation of water temperature higher than 37 degrees C, but it cannot be considered as a typical device for heating the warming of a patient (e.g. patient heating systems). The statement of fulfilment to this standard may be confusing for the user, as in some cases the requirements of this standard are not met. One of such a requirement is the limit of the temperature of applied parts that are in to be in contact with the patients – TrueCool device provides the possibility to achieve 42 degrees of Celsius, and in automatic modes it cannot be changed by the user, where the standard limits the temperature to 40 degrees of Celsius. Due to this requirements several other requirements cannot be achieved – level of triggering of the protection system levels. The standard generally applies to the devices used for the typical heating warming of the patient, but does not fits to the devices where the main function is to keep the patient in hypothermia condition and heating warming is used only to control the level of cooling level (but of course in some aspects, temporarily, the device may behave as typical patient warming system – but still, this may be a short time situation, not like in the typical devices used for patient warming). The blankets/pads/mattresses are not an integral part of the TrueCool device. – they are a separate medical device with its own Technical File. The standard in its cope, clearly indicates, that the standard does not apply to cooling devices, thus, due to the described difficulties and no unacceptable risks, thus, this standard was not included in the documentation, but in many respects its provisions were taken into account (also, included is the temperature stability test performed on applied parts).

EN 60601-2-21:2009	This standard specifies the safety requirements for infant radiant warmers, where TrueCool device is design to use other method of temperature control of the patient's body – conduction instead of radiation, which exist in the typical radiant warmers. Additionally, the description for exclusion of EN 80601-2-35:2009/A1:2016 also applies to this standard.
--------------------	--

18. Transportation and storage of sensors and blankets

Sensors

Store, transport and handle in areas with temperatures above -20 ° C and below 50 ° C.

Blankets

Keep blankets in a dry and clean place.

Storage and transport temperature below 50°C

19. Accessories:

Stop valve

The stop valve is a safety valve that prevents the flow of liquid in the opposite direction.

Working pressure 10 bar

Operating temperature from 0 to + 60 °C

Medium compressed air / Vacuum

Material-polybutylene terephthalate / nickel-plated brass

NBR seals - made of rubber with increased scratch resistance

Under the hose – 12 mm

Length-67 mm

ØP 20.5 mm

Throttle

Valve (needle valve) for adjusting the flow rate of the coolant.

Working pressure -0.95 to + 15 bar

Operating temperature from 0 to + 60 °C

Medium compressed air / Vacuum

Material polybutylene terephthalate / nickel-plated brass


Applicable to PU, PA, PE wires

Thread G1/4

Under the hose 12 mm

L1 49 mm

L2 34 mm

Mono-stereo-test-plug

Plug for checking the selected mode. Single or multiple channel

- 1 pc. Plug; 6.35 mm jack; Dad; Stereo; Straight; on wire; Soldering
- 2 pcs. Plug; 6.35 mm jack; Dad; Mono; Straight; on wire; Soldering
- 3 pcs. Resistor 1kOhm; 0.6 W; $\pm 0.1\%$; 15ppm / $^{\circ}\text{C}$
- 3 pcs. 360 Ohm Resistor; 0.6 W; $\pm 0.1\%$; 15ppm / $^{\circ}\text{C}$

In each connector, there are two resistors connected in a row to get a resistance of 1.36 kOhm, which corresponds to a temperature of about 36.9 $^{\circ}\text{C}$

CST display 7 inch assembly

Color, liquid crystal, touch screen display - touch screen. To display, set parameters, 7"

CST display 10 inch assembly

Color, liquid crystal, touch screen display - touch screen. To display, set parameters, 10"

Bypass valve

The bypass valve is an adjustable pressure relief valve for the cooling / heating fluid.

Damper


A device for absorbing vibration. Dampers are designed for elastic installation, damping of outline, active and passive vibration isolation, noise isolation, force transmission, protective installation, etc. Dampers are made on the basis of general purpose rubber, as well as rubber with increased resistance to scratches, also available with polyester brocade. (SBR/NR/NBR).

D1 - 20 mm

H1 - 20 mm

Screw 1 - M6

Hardness - standard 55 Shore A

Liquid level sensor

The liquid level sensor is required to monitor the level of the cooling / heating liquid in the tank

Installation: Side Entrance

Material: Polypropylene

Side: Outside Side

Feedback form: Form C-single-pole

Contact Material: Ruthenium
Max. switching voltage: 250 VAC or 200V
Maximum Switching Current: 1.0
Maximum Power: 15W
Housing temperature range °C 80
Body Material: Fiberglass Polypropylene
Body Color: Charcoal
Float Material: Fiberglass Polypropylene
Connection type: Stripped cables Cable
type:
2 x 18AWG 32/0.2 mm PVC
insulated 0.5 m long to
UL/CSA/BS6361
Cable Color: Black/Black
Minimum Specific Gravity: 0.7
Operating temperature range °C: 80
Medium / Liquid: Water
Working pressure bar: max. 4 bar
Sealing Components: Silicone

Pump

This is a diaphragm pump, which is a set of flow membranes and a piston mechanism mounted on the engine axis.


W-120 mm D-60 mm
H-130 mm c-48 mm
L-274 mm f-31 mm
a-40 mm g-135 mm
b-64 mm I-4-9x12
c-100 mm
Outlet hole - ¼ inch
Tube -16 mm
Maximum power – 45/52
Standard power – 22/34
Specific gravity-less than 1.1
Number of phases - 1

Schurter cover

Schurter protection / cover is designed to protect the exposed parts of the power cable, as well as the terminals of the connectors of the device from unintentional contact, dust and moisture. It is made of black Polyvinyl chloride and is easy to install.

Compressor

Sealed piston compressor. Creates pressure for the movement of liquid in the tubes

Rated Voltage 220-240 V

Frequency 50 Hz

Condenser

L-250 mm

S-130 mm

LT - 303 mm

M-73 mm

H - 253 mm

ST-203 mm

Weight-4.1 kg

Power - 1650W

Angular

Connecting corner

Option 1

Working pressure -0.95 to 15 bar

Operating temperature from 0°C to +60°C

Medium compressed air / Vacuum

Material polybutylene terephthalate / nickel-plated brass

Applicable to PU, PA, PE wires

Under the hose 12 mm

L1 29.5 mm

L2 13.5 mm

D1 21.5 mm

D2 4.2 mm

Option 2

Gallery-0,99 bar - + 20 bar

Working temperature-20°C to 80°C

Began in a vacuum

Material nickel-Llantium / technopolymer

Material of Colla talm society Nero 304

Profkommerybutelderneftelat 70

Coll cartoon technopolymer

Standard 12mm

A 12.7 mm

B 9,5 mm

L1 31.5 mm

L2 28.5 mm

CH1 16 mm

CH2 25 mm

D 21.5 mm

Option 3

Working pressure -0.95 to 15 bar

Operating temperature from 0°C to +60°C

Medium compressed air / Vacuum

Material PBT plastic / nickel-plated brass

Applicable to PU, PA, PE wires

Thread G1/4

Under the hose 12 mm

L1 7.5 mm

L2 36 mm

L3 29.5 mm

SW 21 mm

D1 21.5 mm

Solenoid body

Mechanical valve with automatic position correction via internal feedback. Designed to adjust the temperature, cooling / heating liquid

Solenoid coil

The solenoid coil is designed to adjust the temperature of the cooling / heating liquid.

Frequency 50Hz

Voltage 24V

Weight-0.5 kg

Steel frame

The steel frame is used to install the hypothermia machine itself

ANNEXES

Annex A: List of recommended of substances for fluid circuit disinfection

The below list contain drug names of which the main active ingredient is Sodium Hypochlorite. If you have concerns about using any of the listed drug, or drug is not available in your area, please contact the manufacturer or use a bleach as an alternative solution.

No.	Name
1	AD Gel
2	Advance 12A (Advance Chemicals)
3	Alcavis Bleach Wipe 1:10 (Alcavis HDC)
4	All-Flex LF (Lalema)
5	All-Flex (Lalema)
6	All Foam (Maxim Chemical)
7	Amuchina (Pisa Laboratorios)
8	Amukline Med (Amuchina)
9	Antiformin
10	B 90 Disinfecting Degrading (Las Produits Agro-B)
11	B814 Detergent (Norchem Industries)
12	Basic 12 (Basic Packaging Industries)
13	Carrel-Dakin
14	Chemspec DFC (Chemical Specialties Manufacturing Corporation)
15	Chlor-O-Fec (Delitek SCS)
16	Chlorasol (Pharmateam Marketing)
17	Chlorasol (SSL International)
18	Chlorine (WestAgro)
19	Chlorobich (Constant America)
20	Chlorox
21	Chlorox
22	Clor-5 (Glen Chemicals)
23	Cloralex
24	Cloropool
25	Clorox Clean-Up Bleach Cleaner (The Clorox Company)
26	Clorox Clean-Up Dilutable Bleach Cleaner (The Clorox Company)
27	Clorox Commercial Solutions Clorox Clean-Up Disinfectant Bleach Cleaner (The Clorox Company)
28	Clorox Commercial Solutions Clorox Clean-Up Disinfectant Dilutable Bleach Cleaner (The Clorox Company)
29	Clorox Commercial Solutions Clorox Professional Disinfecting Bleach Wipes (The Clorox Company)
30	Clorox Commercial Solutions Ultra Clorox Disinfecting Bleach (The Clorox Company)
31	Clorox Disinfecting Toilet Bowl Cleaner with Bleach (The Clorox Company)
32	Clorox Disinfecting Toilet Bowl Cleaner with Bleach-1 (The Clorox Company)
33	Clorox
34	Colix (Prolebor Saude Ocupacional)
35	Comet Cleaner with Bleach (Procter & Gamble)
36	DFC Oxidizing Disinfectant / Disinfectant Cleaner (Cogent Environmental Solutions)
37	DFC Sodium Hypochlorite Disinfectant / Disinfectant Cleaner (Cogent Environmental Solutions)
38	Dakin Cooper Stabilise (Cooperation Pharmaceutique Francaise)
39	Dakin's Solution 0.5% Modified (Century Pharmaceuticals)
40	Dakin's Solution Modified (Century Pharmaceuticals)
41	Dakin's Solution (Century Pharmaceuticals)
42	Dakin (Cooperation Pharmaceutique Francaise)
43	Dakincooper (Melliana)
44	Dakina Solution (Century Pharmaceuticals)
45	Della-San (Delaval)
46	Deosan Green Label Steriliser
47	Desi Forte (Constant America)
48	Desifome II (Constant America)
49	Dilute Sodium Hypochlorite
50	Dispatch
51	Exsept (Alcavis HDC)
52	Fantastik Cream Cleanser with Baking Soda and Disinfecting Bleach (S. C. Johnson and Son)
53	Fantastik for Total Bathroom with Bleach (S. C. Johnson and Son)
54	Fantastik with Bleach All Purpose Cleaner (S. C. Johnson and Son)
55	Formula 12167 (Zep Manufacturing Company)
56	Fresh Scent Clorox Clean-Up Bleach Cleaner (The Clorox Company)
57	G-12 (G.K. Chemical Specialties)
58	Gluco-Chlor (Medrol)
59	Hertel Household Bathroom Disinfectant Cleaner with LA Parisienne Bleach (Lavo)

No.	Name
61	Household Bleach
62	Hydorite
63	Hypercidin (Electramed)
64	Hypure
65	Javel Water
66	Javelle Water
67	Javex 5 Bleach by Clorox (The Clorox Company)
68	Javex 5 Bleach (Colgate-Palmolive)
69	Javex Bleach by Clorox (The Clorox Company)
70	Javex Bleach (Colgate-Palmolive)
71	Javex
72	Klorocin
73	La Parisienne (Lavo)
74	Laundry Super Bleach (Constant America)
75	Lavo 12 (Groupe Lavo)
76	Lavo Pro (Lavo)
77	Liquid Chlorine (West Penetone)
78	Liquido de Dakin (Casa Granado Laboratorio Farmacias e Drogarias)
79	Lysol Brand Disinfectant All Purpose Cleaner with Bleach (Reckitt Benckiser)
80	Lysol Brand Disinfectant Bleach Toilet Bowl Cleaner (Reckitt Benckiser)
81	Lysol Disinfectant All Purpose Cleaner with Bleach (Reckitt Benckiser)
82	Lysol Disinfectant Toilet Bowl Cleaner with Bleach (Reckitt Benckiser)
83	Maranon H (Doloriet)
84	Milton Liquid (Cauta Healthcare)
85	Milton Liquid (Procter & Gamble)
86	Milton (Specialita Igienico Terapeutice)
87	Nacron (Teofarma)
88	Natrum Hypochlorosum (Boliron)
89	Natrum Hypochlorosum (Standard Homeopathic)
90	Neo-Cleaner
91	Neoseptol CL
92	Old Dutch 12 (Groupe Lavo)
93	PCS 1000 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
94	PCS 500 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
95	PCS 5000 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
96	PCS 5000 Oxidizing Disinfectant / Disinfectant Cleaner (Process Cleaning Solutions)
97	PCS Hypochlorite Disinfectant / Disinfectant Cleaner Concentrate (Process Cleaning Solutions)
98	PCS Oxidizing Disinfectant / Disinfectant Cleaner Concentrate (Process Cleaning Solutions)
99	PCS Oxidizing Disinfectant (Process Cleaning Solutions)
100	PCS Sodium Hypochlorite Disinfectant (Process Cleaning Solutions)
101	Parozone
102	Powerfoam C (Asaladin Superior Cleaning Systems)
103	President's Choice Ultra Bleach Disinfectant (KIK Corporation)
104	Princess Auto Powerfoam C (Princess Auto)
105	Purin B
106	Sani-12 (Rochester Midland)
107	Scrubbing Bubbles Bathroom Cleaner (S. C. Johnson and Son)
108	Serca Auto Dish Disinfectant Sanitizer (Sysco)
109	Sodium Hypochlorite (Univar Canada)
110	Soft Scrub with Bleach Cleanser (Dial Canada)
111	Texant
112	Tilex Mold and Mildew Remover (The Clorox Company)
113	Triitol KF (B. Braun)
114	Ultra Clorox Disinfecting Bleach (The Clorox Company)
115	Ultra Javex Professional Bleach by Clorox (The Clorox Company)
116	Virocidin-X (Kam Scientific)
117	X-4 (Ecolab)
118	Youxiolin
119	Zinclin (Westfalla-Surge)

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No.	Наименование
61	Household Bleach
62	Гиклорит
63	Гиперкидин (Electramed)
64	Hypure
65	Javel Water
66	Javelle Water
67	Отбеливатель на основе хлора Javex 5 (компания Clorox)
68	Отбеливатель Javex 5 Bleach (Colgate-Palmolive)
69	Отбеливатель Javex Bleach by Clorox (The Clorox Company)
70	Отбеливатель Javex Bleach (Colgate-Palmolive)
71	Javex
72	Клороцин
73	La Parisienne (Lavo)
74	Laundry Super Bleach (Constant America)
75	Lavo 12 (GroupeLavo)
76	Lavo Pro (Lavo)
77	Liquid Chlorine (West Penetone)
78	Liquido de Dakin (Casa Granado Laboratorio Farmacias e Drogarias)
79	Lysol Brand Disinfectant All Purpose Cleaner with Bleach (Reckitt Benckiser)
80	Lysol Brand Disinfectant Bleach Toilet Bowl Cleaner (Reckitt Benckiser)
81	Lysol Disinfectant All Purpose Cleaner with Bleach (Reckitt Benckiser)
82	Lysol Disinfectant Toilet Bowl Cleaner with Bleach (Reckitt Benckiser)
83	Maranon H (Dolorgiet)
84	Milton Liquid (Ceuta Healthcare)
85	Milton Liquid (Procter & Gamble)
86	Milton (SpecialitalgienicoTerapeutiche)
87	Nacлон (Teofarma)
88	NatrumHypochlorosum (Boiron)
89	NatrumHypochlorosum (Standard Homeopathic)
90	Neo-Cleaner
91	Neosept CL
92	Old Dutch 12 (GroupeLavo)
93	PCS 1000 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
94	PCS 500 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
95	PCS 5000 Oxidizing Disinfectant / Disinfectant Cleaner Wipes (Process Cleaning Solutions)
96	PCS 5000 Oxidizing Disinfectant / Disinfectant Cleaner (Process Cleaning Solutions)
97	PCS Hypochlorite Disinfectant / Disinfectant Cleaner Concentrate (Process Cleaning Solutions)
98	PCS Oxidizing Disinfectant / Disinfectant Cleaner Concentrate (Process Cleaning Solutions)
99	PCS Oxidizing Disinfectant (Process Cleaning Solutions)
100	PCS Sodium Hypochlorite Disinfectant (Process Cleaning Solutions)
101	Parozone
102	Powerfoam C (Aaladin Superior Cleaning Systems)
103	President's Choice Ultra Bleach Disinfectant (KIK Corporation)
104	Princess Auto Powerfoam C (Princess Auto)
105	Purin B
106	Sani-12 (Rochester Midland)
107	Scrubbing Bubbles Bathroom Cleaner (S. C. Johnson and Son)
108	Serca Auto Dish Disinfectant Sanitizer (Sysco)
109	Sodium Hypochlorite (Univar Canada)
110	Soft Scrub with Bleach Cleanser (Dial Canada)
111	Texant
112	Tilex Mold and Mildew Remover (The Clorox Company)
113	Tiutol KF (B. Braun)
114	Ultra Clorox Disinfecting Bleach (The Clorox Company)
115	Ultra Javex Professional Bleach by Clorox (The Clorox Company)
116	Virocidin-X (Kam Scientific)
117	X-4 (Ecolab)
118	Youxiaolin
119	Zincin (Westfalia-Surge)

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