NOTARIAL CERTIFICATE

TO ALL TO WHOM these presents shall come. I, MOHD YACOB KARIM, a Notary Public, officiating at Lot 10, 3rd Floor, Pusat Dagangan UMNO Shah Alam, Persiaran Damai, Seksyen 11, 40100 Shah Alam, Selangor Darul Ehsan, Malaysia duly authorised and appointed, practising in the State of Selangor Darul Ehsan, Malaysia DO HEREBY CERTIFY AND ATTEST THAT I was present on this 24th Day of February, 2017 and did see PN. NOOR AKILAH SAIDIN, the QA DEPUTY GENERAL MANAGER of TOP GLOVE SDN. BHD. (Company No. 220483-T) sign and execute the MEDICAL PRODUCT INFORMATION DATA SHEET dated 24th February 2017 and annexed hereto.

IN FAITH AND TESTIMONY, whereof I have hereunto subscribed my name and affixed my Seal of Office Shah Alam, Selangor Darul Ehsan, Malaysia this 24th Day of February, 2017 in this year of our Lord Two Thousand Seventeen (2017).

MOHD YACOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR DARUL EHSAN
MALAYSIA

MOHD. YACOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR
MALAYSIA
Lot 10, 3rd Floor,
Pusat Dagangan UMNO Shah Alam
Persiaran Damai, Seksyen 11
40100 Shah Alam, Selangor
MALAYSIA
Tel : 603-8513 5610
Fax : 603-8513 5609

24 FEB 2017

APPROVED BY:

QA Deputy General Manager


« TOP GLOVE SDN.BHD. », Malaysia

Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru,

41050 Klang, Selangor D.E., Malaysia

Tel.: + 603 3362 2098 sales@topglove.com.my


Pn. Noor Akilah Saidin,

Medical product information data sheet

Gloves medical diagnostic (examination) disposable non-sterile "BENOVY"

2016

MOHD. YACOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR
MALAYSIA

Lot 10, 3rd Floor,
Pusat Dagangan UMNO Shah Alam
Persiaran Damai, Seksyen 11
40100 Shah Alam, Selangor
MALAYSIA
Tel : 603-5513 5610
Fax : 603-5513 5609

1. Medical product name

Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" have various design variants.

Design variants:

- 1) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 2) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered extra strong (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 3) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powdered (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 4) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination with vanilla flavor (beaded cuff); colour: natural, green, blue, dark blue, pink; size: XS, S, M, L, XL.
- 5) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination (beaded cuff); colour: natural, green, blue; size: XS, S, M, L, XL.
- 6) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating (beaded cuff); colour: natural, green, green mint, blue, black, pink; size: XS, S, M, L, XL.
- 7) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating extra strong (beaded cuff); colour: natural, green, blue, green mint, black; size: XS, S, M, L, XL.
- 8) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination increased strength (long beaded cuff); colour: dark blue; size: XS, S, M, L, XL.
- 9) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with double chlorination (beaded cuff); colour: natural; size: XS, S, M, L, XL.
- 10) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered with green mint flavor (beaded cuff); colour: green mint; size: XS, S, M, L, XL.
- 11) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating with green mint flavor (beaded cuff); colour: green mint; size: XS, S, M, L, XL.
- 12) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with double chlorination increased strength (long beaded cuff); colour: natural; size: XS, S, M, L, XL.
- 13) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with double chlorination (long beaded cuff); colour: natural, blue; size: XS, S, M, L, XL.
- 14) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with single chlorination increased strength (long beaded cuff); colour: blue, light purple, black; size: XS, S, M, L, XL.
- 15) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powder free with polymer coating (beaded cuff); colour: blue, light purple, black, pink; size: XS, S, M, L, XL.

- 16) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powder free with single chlorination (beaded cuff); colour: white, blue, violet blue, light purple, pink, green, red, teal blue, black; size: XS, S, M, L, XL.
- 17) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with single chlorination (beaded cuff); colour: white, blue, orange, black, pink; size: XS, S, M, L, XL.
- 18) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with polymer coating for chemotherapy (beaded cuff/ long beaded cuff); colour: blue; size: XS, S, M, L, XL.
- 19) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powdered (beaded cuff); colour: white, black, blue, violet blue, light purple, green; size: XS, S, M, L, XL.
- 20) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powdered gloves (beaded cuff); colour: white; blue; size: XS, S, M, L, XL.

The Product can be delivered as follows (amount in consumer packaging) :

- 1 pair (2 gloves)
- 5 pairs (10 gloves)
- 10 pairs (20 gloves)
- 20 pairs (40 gloves)
- 25 pairs (50 gloves)
- 50 pairs (100 gloves)
- 100 pairs (200 gloves)

2. Medical product application

2.1.Purpose

Disposable gloves are designed to protect health worker's hands during the examination/treatment of a patient or for other sanitary purposes.

2.2.Scope of use

Gloves are mainly used as a two-sided barrier to protect a patient and health workers from different contaminations.

3. Medical product use

3.1. Indications

Gloves shall be used during diagnostic examinations, therapeutic procedures and medical trials requiring compliance with conditions of patient/health worker protection from cross-infection. When choosing a specific design variant of the gloves it is necessary to consider the following factors to ensure their most comfortable use:

- material the gloves are made of (latex or nitrile),
- cover and coating of external and internal surface (powdered or powder free),
- size,
- utility and texture (non-textured, textured on fingers, textured)
- comfort (roller cuff, cuff without a roller, extended cuff).
- safety: may cause allergic reactions (latex, powdered) or may not cause allergic reactions

(nitrile, powder free).

Attention! The gloves do not provide 100% hand protection from perforation and cuts. The gloves are designed for single use only.

3.2. Contraindications

- - latex gloves may cause allergic reaction
- - powdered gloves may cause allergic reaction


3.3. Medical product application rules

- Every health worker aware of all features of his/her work should sum up individual requirements for the gloves taking into account:
 - actual occupational hazards (evaluation of biological, mechanical and chemical risks);
 - sanitary and epidemiological conditions (sanitary demand);
 - operational features (types of performed procedures and actions, physical efforts, duration of work in gloves; tactile sensitivity need, etc.);
 - skin health of the health worker, allergological history of the health worker and a patient and some other factors.
- Gloves are no substitute for proper hand hygiene. Handwash should be performed in accordance with the mode of procedures.
- It is necessary to choose gloves of appropriate size.
- Gloves should be chosen depending on procedure type, they should be worn easily, should not hamper movements or compress hands.
- When putting on gloves a person should not make great efforts.
- Before putting on gloves for medical procedures, working area should be prepared: sterile instruments and consumables should be put on a manipulation table ready for use. It is necessary to ensure their accessibility without breaching the sterility of hands.
- Using additional pair of gloves may reduce gloves damage rate when working with sharp instruments by 2-4 times.
- After gloves removal it is necessary to wash hands with soap and disinfect them with dermal antiseptic.

3.4. Putting on gloves

It is recommended to follow the techniques described below in order to avoid external surface contamination when putting on the gloves. There are 2 techniques to put on the gloves: open and closed. Whatever technique is chosen, the most important is to avoid external surface contamination.

Open technique steps

- Open package and take two gloves with your finger tips of the left hand from inside of the left glove, so that the thumb tip is inside the right glove, and hold them together vertically. Make sure the fingers contact only with an internal surface of the gloves.
- Pass the right hand through the cuff holding the gloves in your left hand.
- Pull the folded glove until your right hand is inside the glove.
- Put right hand fingers into the left glove cuff and draw it on the left hand.
- Draw the left hand glove cuff on a medical coat cuff with four fingers.
- Put your left hand fingers in the glove under a folded glove cuff and draw it on a medical coat cuff.
- Draw on the cuff with all four fingers and fix the gloves for comfort use.
- In case of powered gloves, the powder should be aseptically removed from the glove surface.

Close technique steps

- Open secondary package before washing your hands.

- Draw coat cuffs on your hands, grab the left glove with the right hand. Hold the hand supine, fingers straight. Put the glove on the left wrist so that its fingers are pointed at you and the thumb lays on the wrist. Press down the cuff by your left thumb.
- Put the thumb on the end inside the cuff. Jam your right hand and draw the glove on the fingers.
- Straightening out your left hand fingers and holding the cuff in the right hand, draw on the glove with four right hand fingers till it is drawn on the left hand fingers.
- Repeat the steps with the second glove. In particular: Put the right glove on the right wrist using the left hand in the glove so that its fingers are pointed at you and the thumb lays on the wrist. Press down the cuff by the right thumb.
- Put the left thumb on the end inside the cuff. Jam your right hand and draw the glove on the fingers.
- Straightening your right hand fingers and covering the wrist by the cuff of the medical coat draw the glove on with four left hand fingers till it is drawn on the right hand fingers.
- Now, when the gloves are put on both hands they can be fixed without the risk of contamination.
- To avoid gloves slipping from the coat cuffs, please ensure the sufficient length of the medical coat sleeves and sufficient cuff length up to the heel of the hand, as well as avoid folds on the glove cuff.
- In case of powdered gloves, the powder should be aseptically removed from the glove surface.

3.5. Requirements for working in the gloves

- The gloves should be changed when switching from “dirty procedures” to “clean procedures” or going from one patient to another.
- It is necessary to change the gloves:
 - if any visual defect is observed;
 - after potential or actual perforation by a needle, suture material, bone, etc.;
 - after accidental electric shock to a hand from an electrosurgical unit;
 - if a glove starts to bloat, strain, stick or stays loosely on the hand due to absorption of liquids or lubricants;
- In case of using powdered gloves it is necessary to aseptically remove powdering agent by: washing it off with sterile water or washing hands in a sterile tank or using a sterile drape. Powder free gloves do not require additional treatment.
- The gloves can not be cleaned during procedures. They should be changed if necessary.
- Damaged or perfused gloves can not be used. If a glove is damaged, both gloves should be changed immediately because it is impossible to remove one glove without contamination of another.
- Remove the gloves carefully avoiding hand skin contamination and spillage of liquids.
- It is not recommended to clean disposable latex gloves with alcohol-containing solutions, because these solutions may damage a protective cover of the gloves. Gloves made of synthetic materials (nitrile, etc.) are proof against the action of antiseptic chemical agents
- Strictly observe safety precautions while working with sharp piercing objects and instruments.

3.6. Removing the gloves

It is recommended to comply with the following glove removal algorithm:

- Take the glove near the heel of the hand with the fingers of the other hand.
- Remove the glove and hold it with the other hand.
- Then put the fingers of your free hand into the glove cuff on another hand and remove the glove as the first one.
- Turn the glove inside out and use it as a bag for the second glove.
- Put the gloves in a special container. Then clean the hands.

4. Medical product description and its specification.

4.1. Depending on the used material two types of gloves are manufactured:

- type 1 – natural rubber latex gloves;
- type 2 – synthetic rubber (nitrile) gloves

4.2. There are the following types of gloves depending on the cover:

- texture of some part or the whole glove surface;
- non-textured surface;
- powdered surface;
- non-powdered surface;
- single chlorination;
- double chlorination;
- with inner polymer coating.

4.3. Flexicoat 88 polymer is used for internal coating of gloves. This polymer coating ensures additional durability of gloves, isolates skin from natural latex proteins and considerably facilitates the putting on, as well as prevents sticking of gloves.

Gloves can be processed using gaseous chlorine (chlorination). The chlorination can be single or double (both sides).

USPXXI grade modified cornstarch is used for gloves powdering.

4.4. Gloves have five nominal sizes:

- extra small (XS), size code 6 or below;
- small (S), size code 6.5;
- medium (M), size code 7, 7.5;
- large (L), size code 8, 8.5;
- extra large (XL), size code 9 and above.

4.5. The following dyes are used during gloves manufacture:

Green Mint Colour – Colortrend Green 1211/79

Blue colour – L7100

Black colour – Farsperse Black

Green colour- Farsperse Green PG 708R

Pink colour – Esperse L4600
Dark blue colour – Esperse Blue L7100
Light purple colour - Esperse L8100C
White colour – NIL
Orange colour - Farsperse Orange MP 01411
Teal blue colour – MPL 1244
Red colour - Pantone 7635 C
Violet blue colour – pantone 2726 C

4.6. The following aromas are used:
Mint flavour. 50% Peppermint Oil
Vanilla flavour: 13503, Matrix Oleochem

4.7. High-impact (increased strength) gloves feature increased side thickness, textured surface and higher resistance to damage.

High-impact (increased strength) gloves are used in conditions of a high infection and contamination risk during procedures associated with the possibility of mechanical damage and requiring reliable hands protection.

4.8. Chemotherapy gloves are resistant to permeation of chemotherapy (cytotoxic) drugs. The gloves are tested for protection against cytotoxic drugs according to standard ASTM D 6978-05: Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs

The gloves are designed to prevent infection and protect hands of health care professionals against chemotherapy drugs.

4.9. Extra-strong gloves are made from latex and can be powder free with inner polymer coating and powdered, this gloves differ from conventional similar type of gloves in performance glove thickness wall. For smooth surface at least 0.10 mm, for textured surface at least 0,13 mm.

4.10. Gloves with green mint aroma have a light green colour and a discreet odour of mint.

4.11. Gloves with vanilla aroma have a natural colour and a discreet odour of vanilla.

4.12. The products are designed for single use only.

4.13. Class according to the degree of the potential usage risk –1.

5. Specification

Width, length and thickness of gloves, depending on the size, shall correspond to values specified in Table 1.

5.1.If measured according to the points (shown in Figure 1), glove size should comply with width and length (see Table 1) in case of inspection level and acceptance quality level (AQL) specified in Table 2.


Figure 1 - Points to measure width and length of the glove

Table №1 Gloves sizes

Size	Palm's width, mm	Minimum length, mm	Minimum length (extended roller cuff), mm	Minimum thickness at measurement points (Figure 2), mm
6.0 or less/ Extra small (X-S)	$\leq 82/\leq 80$	220	300 ± 10	0.08 — for smooth sections 0.11 — for textured sections
6.5/ Small (S)	$83\pm 5/80\pm 10$	220	300 ± 10	
7.0/ Medium (M)	$89\pm 5/95\pm 10$	230	300 ± 10	
7.5 /Medium (M)	$95\pm 5/95\pm 10$	230	300 ± 10	
8.0 /Large (L)	$102\pm 6/110\pm 10$	230	300 ± 10	
8.5 /Large (L)	$109\pm 6/110\pm 10$	230	300 ± 10	
9.0 or more/ Extra large (X-L)	$\geq 110/\geq 110$	230	300 ± 10	

The length is measured as the shortest distance from the tip of the middle finger to the edge of the cuff.

The width is measured on a flat surface using the midpoint between the base of the index finger and the thumb.

5.2. Thickness of two sides of an intact glove is measured using a thickness gage with pressure on a base of (22 ± 5) kPa in three points shown in the Figure 2 at a distance of (13 ± 3) mm from the tip of the middle finger, in the centre of a palm, and at a distance of (25 ± 5) mm from the edge of the cuff. Thickness of one side in each point shall be equal to a half of thickness value for two sides and has to meet the standard given in Table 1 in case of inspection level and acceptance quality level (AQL) specified in Table 2.

If, in case of visual inspection, thinning-out in any section of a glove is revealed, thickness of one side of a glove in this section is measured; it shall be at least 0.10 mm for non-textured and at least 0.13 mm for textured surface.

Thickness of glove edge shall be not more than 2.50 mm.


Figure 2

1-finger tip, 2-palm, 3-cuff edge

A= (13 ± 3) mm, B= (33 ± 5) mm, C= (48 ± 9) mm

High-impact (increased strength) gloves Cuff length: 310 ± 10 mm.

Palm section thickness (centre of a palm): $0.20 \text{ mm} \pm 0,02 \text{ mm}$

Finger section thickness (13 mm from a tip): $0.23 \text{ mm} \pm 0,02 \text{ mm}$

For the purpose of control, gloves are selected and checked in accordance with ISO 2859 -1. The inspection levels and acceptance quality levels (AQL) for gloves are provided in Table 2.

2 Inspection levels and acceptance quality levels (AQL)

Property	Test defect	Inspection level	AQL
Size (length, thickness, width)	Width, length, thickness	S-2	4.0
Hermeticity	Presence of holes	G-1	2.5
Force at break and elongation at break (before and after accelerated ageing):	Before and after accelerated ageing	S-2	4.0
Physical and mechanical parameters	Before/after accelerated ageing	S-2	4.0
Residual powdering agent	Exceeds the recommended limit	N=5	N/A
Powdering agent content	Exceeds the recommended limit	N=2	N/A
Water soluble protein content (for latex only)	Exceeds the recommended limit	N=3	N/A

6. Storage and transportation conditions:

6.1 The product may be transported by any means of transport, in accordance with transport regulations in force for this type of transport, in its original transport packages (pressboard boxes) protected from direct sunlight.

Gloves shall retain their properties under operating conditions in case of the following environmental condition: temperature and relative humidity of air:

temperature – +10-+35 °C;

relative air humidity – up to 80% at +25 °C.

Gloves, packed into transport containers, shall be resistant to the following climatic effects:

temperature – -30-+50 °C;

relative air humidity – up to 98% at 35 °C, but without moisture condensing.

6.2. The gloves shall be stored in original multiple pressboard packages in a well-ventilated room away from heat sources, direct sunlight, direct exposure of light and electrical equipment. Condensation may occur if storage temperature is less than 10 °C. If storage temperature is less than 10 °C, the gloves shall be warmed up to 25-30 °C before use, otherwise difficulties may occur due to elasticity degradation of the gloves.

6.3. Store in a dry place protected from direct sunlight and high humidity at 10 °C-30 °C.

6.4. Store the gloves away from sources of ignition.

7. Packaging

7.1. Packaging shall provide protection from mechanical and climatic factors during transportation and storage, as well as ensure complete vehicle load (capacity) use and convenience of handling operations.

7.2. Packing of products shall conform to the requirements of the manufacturer and terms of the contract between the manufacturer and the authorized representative.

Consumer packaging

Packages with gloves are packed in cardboard boxes. Consumer packaging (outer package) protects gloves from external mechanical damage and allows to visually determine the size of the gloves.

Transport package

Consumer packages with gloves are packed in corrugated cardboard boxes. Cardboard boxes shall be durable enough to ensure the quality of the product during transportation and storage.

Instructions for use are printed on the consumer packaging.

8. Product set.

The following product sets are offered:

Amount in consumer packaging:	Amount in transport packaging
• 1 pair (2 gloves)	500 pairs (1000 gloves)
• 5 pairs (10 gloves)	500 pairs (1000 gloves)
• 10 pairs (20 gloves)	500 pairs (1000 gloves)
• 20 pairs (40 gloves)	500 pairs (1000 gloves)
• 25 pairs (50 gloves)	250 pairs (500 gloves) / 500 pairs (1000 gloves)
• 50 pairs (100 gloves)	500 pairs (1000 gloves)
• 100 pairs (200 gloves)	1000 pairs (2000 gloves)

Number of pairs in transport packages may vary in case it is agreed upon with the authorized representative and is in accordance with the requirements of the contract concluded between manufacturer and authorized representative.

9. Medical product supply package.

- Medical examination disposable non sterile gloves "BENOVY" made in one of the design variants.

Design variants:

- 1) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 2) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered extra strong (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 3) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powdered (beaded cuff); colour: natural, green, blue, green mint, dark blue, pink; size: XS, S, M, L, XL.
- 4) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination with vanilla flavor (beaded cuff); colour: natural, green, blue, dark blue, pink; size: XS, S, M, L, XL.
- 5) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination (beaded cuff); colour: natural, green, blue; size: XS, S, M, L, XL.
- 6) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating (beaded cuff); colour: natural, green, green mint, blue, black, pink; size: XS, S, M, L, XL.
- 7) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating extra strong (beaded cuff); colour: natural, green, blue, green mint, black; size: XS, S, M, L, XL.
- 8) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with single chlorination increased strength (long beaded cuff); colour: dark blue; size: XS, S, M, L, XL.

- 9) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with double chlorination (beaded cuff); colour: natural; size: XS, S, M, L, XL.
- 10) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex smooth powdered with green mint flavor (beaded cuff); colour: green mint; size: XS, S, M, L, XL.
- 11) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with inner polymer coating with green mint flavor (beaded cuff); colour: green mint; size: XS, S, M, L, XL.
- 12) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" from natural latex textured powder free with double chlorination increased strength (long beaded cuff); colour: natural; size: XS, S, M, L, XL.
- 13) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with double chlorination (long beaded cuff); colour: natural, blue; size: XS, S, M, L, XL.
- 14) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with single chlorination increased strength (long beaded cuff); colour: blue, light purple, black; size: XS, S, M, L, XL.
- 15) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powder free with polymer coating (beaded cuff); colour: blue, light purple, black, pink; size: XS, S, M, L, XL.
- 16) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powder free with single chlorination (beaded cuff); colour: white, blue, violet blue, light purple, pink, green, red, teal blue, black; size: XS, S, M, L, XL.
- 17) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with single chlorination (beaded cuff); colour: white, blue, orange, black, pink; size: XS, S, M, L, XL.
- 18) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powder free with polymer coating for chemotherapy (beaded cuff/ long beaded cuff); colour: blue; size: XS, S, M, L, XL.
- 19) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile finger textured powdered (beaded cuff); colour: white, black, blue, violet blue, light purple, green; size: XS, S, M, L, XL.
- 20) Gloves medical diagnostic (examination) disposable non-sterile "BENOVY" nitrile textured powdered gloves (beaded cuff); colour: white; blue; size: XS, S, M, L, XL.

- 1 information data sheet (or instructions indicated on the package)

- 1 package (1 pair (2 gloves)/ 5 pairs (10 gloves) /10 pairs (20 gloves)/ 20 pairs (40 gloves)/ 25 pairs (50 gloves)/ 50 pairs (100 gloves)/ 100 pairs (200 gloves))

10. Safety requirements

- Materials the gloves are made of shall be toxicologically safe.
- In case of powered gloves, the powder should be aseptically removed from the glove surface.
- The gloves are designed for single individual use only.
- The gloves shall not guarantee 100 % protection from punctures and cuts.

11. Side effects

Gloves made of natural rubber latex can cause allergic reactions.

Powdered gloves can cause allergic reactions.

Health workers can notice the following negative effects caused by the use of gloves: non-allergic contact dermatitis, allergic contact dermatitis and immediate hypersensitivity reactions presented by urticaria, pruritus and systematic allergic reactions.

Allergic reactions are generally caused by allergy to latex, i.e. personal sensitivity reactions of health workers to chemical components of latex as latex is the most popular material used to produce medical gloves. Use of powder to facilitate the process of putting on gloves also raises the risk of health workers' allergy and can cause some postoperative complications experienced by patients.

12. Medical product application conditions

Products shall be used in healthcare institutions, means of transport as well as in the field and at home in situations when hands of health workers need protection during examination/treatment of patients as well as during other medical procedures. Gloves are mainly used as a double-sided barrier to protect patients and healthcare workers from various contaminations.

The gloves can be used as standard safety measures.

No specific qualification or education shall be required to use the gloves.

The following application conditions shall be observed:

Personal use.

Single use (disposable product).

Kind of contact with human organism: short-time contact with undamaged skin and mucous membrane, biological liquids.

Required room temperature and relative humidity:

temperature: +10 - +35 °C; relative humidity: up to 80% if the temperature is +25 °C.

Apply antiseptic medicine on your hands before putting on gloves depending on planned operations and procedures and remove powder agent maintaining sterility (for powdered gloves).

Put on a glove accurately on your hand, carefully expand the cuff up to your wrist, expand accurately the glove on your hand with the other hand. To take off gloves make a fold on the left glove with fingers (use your right hand) touching only external part and take off your left-hand glove turning it inside out and holding the fold of the glove. With the left-hand glove in your right hand hold the right-hand glove by the fold (on the inside), take off the right-hand glove turning it inside out: the left-hand glove is in the middle of the right-hand glove. Put the gloves in the dustbin. Do not reuse them.

When the gloves are taken off wash your hands and apply antiseptic!

Do not reuse gloves. Dispose of gloves in accordance with the specified procedure for disposal of medical products.

13. Labelling.

Packaging labelling shall conform to the requirements of this document and applicable legislation. The information on packaging is printed in Russian and, if necessary, may be repeated in other languages.

Consumer package labelling shall include:

1. Name of medical product.
2. Name of the manufacturer and manufacturer's legal representative in the Russian Federation.
3. Address of the manufacturer and manufacturer's legal representative in the Russian Federation.

4. Trademark of the manufacturer and/or manufacturer's legal representative in the Russian Federation.
5. Indication of material used.
6. Product colours.
7. Fragrance (for products with fragrance).
8. Chlorination (for products with single or double chlorination).
9. Polymeric cover (for products with polymeric cover).
10. Words "Textured" or "Non-textured", "Powdered" or "Powder free".
11. Size information.
12. Batch number.
13. Number of items in a package.
14. "Date of manufacture" as well as year and month of manufacture.
15. Words "Examination gloves", Words "Disposable".
16. Words "The product is made of natural rubber latex, it may cause allergic reactions" (for products made of natural latex); "Powdered gloves may cause allergic reaction" (for powdered gloves).
17. In case of using powdering agents it is necessary to specify that before using the gloves the powdering agent shall be aseptically removed.
18. Information on cuffs.
19. Shelf life.
20. Storage conditions.
21. Part number.
22. Bar code.
23. Conformity assessment information.
24. Indication of regulatory document in force in a country of the legal representative (Russian Federation).
25. Number and date of marketing authorisation necessary for product distribution in a country of the legal representative (Russian Federation).
26. Medical product designation.
27. Information needed for use, storage and disposal of the product.
28. It is possible to use the symbols for medical products of ISO 15223-1:2012 "Medical devices - Symbols to be used with medical device labels, labelling and information to be supplied - Part 1: General requirements" (and other symbols approved in the country of medical product distribution) (if necessary).

14. Warranty and shelf life

The manufacturer guarantees integrity of all parameters and characteristics of the gloves throughout the shelf life (5 years), if packages are not damaged and storage, transportation and operating instructions are followed.

15. Disposal and destruction procedures

The system of collection, temporary storage and transportation of medical waste should consist of the following steps:

- internal waste collection in a healthcare and/or pharmaceutical organization;
- transportation of waste from the divisions and its temporal storage on the territory of the organization generating waste;

- disinfection/neutralization;
- transportation of waste from the territory of the organization generating waste;
- dumping or disposal of medical waste.

Chemical disinfection waste at place of its generation is performed as a mandatory interim measure in the absence of medical waste handling area of healthcare and/or pharmaceutical organizations, or in the absence of a centralized system of medical waste neutralization in the area.

Methods of safe disinfection and/or neutralization waste shall be chosen depending on the capacity and profile of the healthcare organization, availability of facilities for disinfection/neutralization of waste, methods of waste disinfection/neutralization adopted in the area (incineration, removal, recycling).

Collection, temporal storage and removal of waste shall be performed in accordance with a medical waste handling plan of a healthcare and/or pharmaceutical organization.

16. Medical product manufacturer data

Manufacturer:

«TOP GLOVE SDN. BHD.», Malaysia

Address: Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor D.E.,
Malaysia

Manufacturer and production site address:

«TOP GLOVE SDN. BHD. », Malaysia

Address: Lot 4969, Jalan Teratai, Batu 6, Off Jalan Meru, 41050 Klang, Selangor D.E.,
Malaysia

17. Claims

For all issues related to gloves medical diagnostic (examination) disposable non-sterile "BENOVY" manufactured by TOP GLOVE SDN. BHD, Malaysia, please contact the legal representative of the manufacturer in the Russian Federation: LLC «BENOVY», legal and postal address: lit. J, 1, 5 Verkhniy pereulok, Saint Petersburg, 194292 Tel.: (812) 309-16-32. E-mail: benovy@benovy.ru.


MOHD. YUSOB KARIM
NOTARI AWAM
NOTARY PUBLIC
SELANGOR
MALAYSIA

Lot 10, 3rd Floor,
Pusat Dagangan UMNO Shah Alam
Persiaran Damai, Seksyen 11
40100 Shah Alam, Selangor
MALAYSIA
Tel : 603-5513 5610
Fax : 603-5513 5609


24 FEB 2017

1. Medical product information data sheet
(Total: 16 pages)

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