
Medtronic, Inc.
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
tel 001-763-514-4000
fax 001-763-514-4879

Declaration Letter

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Acheve, Arctic Front, Arctic Front Advance, Arctic Front Advance Pro, CryoConsole, FlexCath, Medtronic

Explanation of symbols • Vysvětlení symbolů • Symbolforklaring • Erläuterung der Symbole • Ερμηνεία των συμβόλων • Explicación de los símbolos • Explication des symboles • A szimbólumok magyarázata • Spiegazione dei simboli • Verklaring van de symbolen • Symbolforklaring • Objasnienie symboli • Explicação dos símbolos • Explicação dos símbolos • Объяснение символов • Vysvětlivky k symbolom • Forklaring av symboler • Sembollerin açıklaması

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1 Description


The Arctic Front Advance Pro Cardiac Cryoablation Catheter (the catheter or the Arctic Front Advance Pro Cryoballoon) is a flexible, over-the-wire balloon catheter used to ablate cardiac tissue. It is used together with a compatible Medtronic 12 Fr inner diameter sheath (the sheath), the CryoConsole, and related components. For device compatibility questions, contact Medtronic Technical Support.

The balloon reaches cryoablation temperatures when refrigerant is injected from the CryoConsole to the balloon segment. A thermocouple positioned inside the balloon provides temperature reading capability. The catheter is introduced into the vasculature by traditional, minimally invasive techniques.

There are two radiopaque markers on the catheter to confirm the position of the balloon using fluoroscopy. The proximal radiopaque marker is located approximately 10 mm (0.394 in) proximal to the balloon. The distal radiopaque marker is located at the end of the injection tube.

Note: The 12 Fr FlexCath sheath is compatible with the catheter. There is one radiopaque marker on the FlexCath sheath, located approximately 5 mm (0.197 in) proximal to the tip of the sheath (see Figure 1). When the catheter's proximal radiopaque marker and the FlexCath sheath's radiopaque marker are aligned, the balloon is approximately 5 mm (0.197 in) outside of the FlexCath sheath. There are two shaft markers on the proximal section of the catheter shaft to visually confirm the position of the balloon within the FlexCath sheath. When the distal marker on the shaft of the catheter is aligned with the handle of the FlexCath sheath, the balloon segment is located inside the FlexCath sheath. When the proximal marker on the shaft of the catheter is aligned with the handle of the FlexCath sheath, the balloon segment is outside the FlexCath sheath (see Figure 1).

Figure 1. Arctic Front Advance Pro Cardiac Cryoablation Catheter

- | | |
|---------------------------------------|----------------------------------|
| 1 Distal catheter radiopaque marker | 5 Proximal catheter shaft marker |
| 2 Deflated balloon segment | 6 Distal catheter shaft marker |
| 3 Proximal catheter radiopaque marker | 7 Deflection mechanism |
| 4 FlexCath sheath radiopaque marker | 8 Blue push button |

The Arctic Front Advance Pro Cardiac Cryoablation Catheter is available in 2 models, as described in the following table:

Model	Inflated balloon diameter
AFAPRO23	23 mm (0.91 in)
AFAPRO28	28 mm (1.10 in)

For details about the CryoConsole and how to use it with the catheter to perform cryoablation procedures, see the *CryoConsole Operator's Manual*.

1.1 Contents of package

The catheter is supplied sterile. The package contains the following items:

- 1 Arctic Front Advance Pro Cardiac Cryoablation Catheter
- product documentation

2 Indications for use

The Arctic Front Advance Pro Cardiac Cryoablation Catheter is indicated for the treatment of patients with atrial fibrillation.

3 Contraindications

The Arctic Front Advance Pro Cardiac Cryoablation Catheter is contraindicated as follows:

- in the ventricle because of the danger of catheter entrapment in the chordae tendineae
- in patients with active systemic infections
- in conditions where the manipulation of the catheter within the heart would be unsafe (for example, intracardiac mural thrombus)
- in patients with cryoglobulinemia
- in patients with one or more pulmonary vein stents

4 Warnings and precautions

Anticoagulation therapy – Administer appropriate levels of peri-procedural anticoagulation therapy for patients undergoing left-sided and transseptal cardiac procedures. Administer anticoagulation therapy during and post-procedure according to the institution's standards. The Arctic Front Advance Pro Cryoballoon was not studied for the safety of changes in anticoagulation therapy in patients with atrial fibrillation.

Balloon inflation or deflation – Do not inflate the balloon inside the sheath. Always verify with fluoroscopy or other appropriate visualization techniques that the balloon is fully outside the sheath before inflation to avoid catheter damage.

- Do not inflate the balloon while the catheter is positioned inside a pulmonary vein. Always inflate the balloon in the atrium and then position it at the pulmonary vein ostium; inflating the balloon in the pulmonary vein may result in vascular injury.
- If the balloon cannot be inflated or deflated using the CryoConsole, have a Manual Retraction Kit on hand during the procedure. (Refer to the *CryoConsole Operator's Manual* for more detailed instructions on the Manual Retraction Kit).

Biohazard disposal – Discard all used catheters and sterile components in accordance with hospital procedures.

Cardioversion or defibrillation during ablation procedure – Disconnect the catheter's electrical connection before cardioversion or defibrillation. Failure to do so may trigger system messages indicating a need for catheter exchange.

Catheter handling – Use extreme care when manipulating the catheter. Lack of careful attention may result in injury such as perforation or laceration.

- Do not use excessive force to advance, withdraw, or rotate the catheter, especially if resistance is encountered. Excessive force may lead to catheter damage, including kinking of the guide wire lumen within the balloon segment.
- Do not use the catheter if it is kinked, damaged, or cannot be straightened.
- Straighten the shaft before inserting or withdrawing the catheter.
- Do not at any time preshape or bend the catheter shaft or balloon segment. Bending or kinking the catheter shaft may damage internal structures and increase the risk of catheter failure. Prebending of the distal curve may damage the catheter.
- Catheter advancement should be performed using fluoroscopy or other appropriate techniques.
- Do not position the cryoballoon catheter within the tubular portion of the pulmonary vein to minimize phrenic nerve injury and pulmonary vein stenosis.

Catheter integrity – Do not use the catheter if it is kinked or damaged. If the catheter becomes kinked or damaged while in the patient, remove it and use a new catheter. Before injecting, the physician should ensure that there is no kink in the catheter.

Circular mapping catheter compatibility – Use only Medtronic circular mapping catheters compatible with the inner lumen of the Arctic Front Advance Pro Cryoballoon. Use of another mapping catheter may damage the catheter or compromise the procedure.

Contrast media – Use appropriate levels of contrast media in patients with comorbidities such as recent history of renal disease. Follow contrast labeling and institutional procedures regarding the appropriate medical strategies to minimize risk when using contrast media.

Correct guide wire or circular mapping catheter insertion and positioning – Do not advance the balloon beyond the guide wire or circular mapping catheter to reduce the risk of tissue damage.

- Ensure that the guide wire or circular mapping catheter is inserted into the catheter and through the balloon portion for adequate support during vascular access insertion. Failure to do so may result in catheter damage.

Cryoablation near prosthetic heart valves – Do not pass the catheter through a prosthetic heart valve (mechanical or tissue). The catheter may become trapped in the valve, damaging the valve and causing valvular insufficiency or premature failure of the prosthetic valve.

Cryoablation – Do not pull on the balloon catheter, circular mapping catheter, sheath, umbilical cables, or CryoConsole while the balloon catheter or circular mapping catheter are frozen to tissue. This may lead to tissue injury. Before moving these components, use appropriate techniques to ensure that the balloon catheter and circular mapping catheter are not adhered to tissue.

Damage to lung or tracheobronchial tree – Damage to the lung or tracheobronchial tree has been observed in some subjects who have undergone left atrial ablation with the Arctic Front family. The physician should consider appropriate medical strategies to minimize the risk of damage to the lung or tracheobronchial tree.

Do not resterilize – Do not resterilize this device for the purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.

Embolism risk – Introducing any catheter into the circulatory system entails the risk of air or gas embolism, which may occlude vessels and lead to tissue infarction with serious consequences. Always advance and withdraw components slowly to minimize the vacuum created and therefore minimize the risk of air embolism.

Environmental limits – Perform cryoablation procedures only within the environmental parameters. Operating outside these parameters may prevent the start or completion of a cryoablation procedure. Refer to Chapter 7, "Specifications", page 8 for environmental parameters.

Esophageal injury – Esophageal ulcerations have been observed in some subjects who have undergone left atrial ablation with the Arctic Front family. As with other forms of left atrial ablation, the physician should consider appropriate medical strategies to minimize the risk of esophageal injury.

Fluid incursion – Do not expose the catheter handle or coaxial and electrical connectors to fluids or solvents. If these components get wet, the system may not function properly. This may lead to patient injury.

Fluoroscopy required for catheter placement – The use of fluoroscopy during catheter ablation procedures presents the potential for significant x-ray exposure to both patients and laboratory staff. Extensive exposure may result in acute radiation injury and increased risk for somatic and genetic effects. Only perform catheter ablation after giving adequate attention to the potential radiation exposure associated with the procedure, and taking steps to minimize this exposure. Give careful consideration before using the device in pregnant women.

For single use only – This device is intended only to be used once for a single patient. Do not reprocess or resterilize this device for the purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Frequent flushing of the guide wire lumen – Flush the guide wire lumen before initial insertion and then frequently throughout the procedure to prevent coagulation of blood in the lumen. Flush the guide wire lumen with heparinized saline after each contrast injection.

Guide wire compatibility – Use only 0.081 cm (0.032 in) or 0.089 cm (0.035 in) guide wires with the catheter. Using another guide wire may damage the catheter or compromise the procedure.

Improper connection – Do not connect the cryoablation catheter to a radiofrequency (RF) generator or use it to deliver RF energy. This may cause catheter malfunction or patient injury.

Induced arrhythmias – Catheter procedures may mechanically induce arrhythmias.

Leakage current from connected devices – Use only isolated equipment (IEC 60601-1 Type CF equipment, or equivalent) with the CryoConsole and catheters or patient injury or death may occur.

Other catheters, devices, or wires – Avoid catheter entanglement with other catheters, devices, or wires. Such entanglement may necessitate surgical intervention.

Phrenic nerve impairment – To reduce the potential for phrenic nerve impairment, perform the following steps:

- Position the balloon as antral as possible and not in the tubular portion of the pulmonary vein. To ensure proper catheter position, observe the balloon shape and the balloon position within the anatomy.
- Monitor the status of the phrenic nerve continuously during right-sided pulmonary vein applications using an appropriate monitoring technique. One common monitoring technique is to continuously pace the phrenic nerve throughout each cryoablation application of the right pulmonary veins. While pacing, monitor diaphragm contraction by placing a hand on the abdomen to assess for loss of capture or changes in the strength of the diaphragmatic contraction.
- Stop ablation immediately if phrenic nerve impairment is observed.

Post-ablation period – Closely monitor patients undergoing cardiac ablation procedures during the post-ablation period for clinical adverse events.

Pressurized refrigerant – The catheter contains pressurized refrigerant during operation. Release of this gas into the circulatory system due to equipment failure or misuse may result in gas embolism.

Pulmonary vein narrowing or stenosis – Catheter ablation procedures inside or near pulmonary veins may induce pulmonary vein narrowing or stenosis. Do not ablate in the tubular portion of the pulmonary vein. The occurrence of this complication may necessitate percutaneous angioplasty or surgical intervention.

Qualified users – This cryoablation system should be used only by or under the supervision of physicians trained in cryoablation procedures.

Required use environment – Cryoablation procedures should be performed only in a fully equipped facility.

RF ablation – Before powering up an RF generator or applying RF energy, disconnect the cryoablation catheter from the CryoConsole to avoid an error message and unnecessary catheter replacement.

Septal damage – Always deflate the balloon and withdraw the balloon into the transseptal sheath before removing the balloon from the left atrium. Crossing the septum while the balloon is unsheathed, inflated, or inflating in the septal puncture site may cause serious septal damage.

Steerable sheath compatibility – Use only a compatible Medtronic 12 Fr inner diameter sheath with the Arctic Front Advance Pro Cryoballoon. Using another sheath may damage the catheter or balloon segment.

Sterile package inspection – Inspect the sterile packaging and catheter before use. If the sterile packaging or catheter is damaged, do not use the catheter. Contact your Medtronic representative.

System compatibility – Use only Medtronic cryoablation catheters, refrigerant tanks, and components with the CryoConsole. The safety and use of other catheters or components has not been tested.

5 Adverse events

Potential adverse events associated with cardiac catheter cryoablation procedures include, but are not limited to, the following conditions:

- Access site complications (e.g. bruising, ecchymosis)
- Anemia
- Anxiety
- Arrhythmia (e.g. atrial flutter, bradycardia, heart block, tachycardia)
- Back pain
- Bleeding from puncture sites
- Bronchial constriction
- Bronchial fistula
- Bronchitis
- Bruising
- Cardiac tamponade
- Cardiopulmonary arrest
- Cerebral vascular accident
- Chest discomfort/pain/pressure
- Cold feeling
- Coronary artery spasm
- Cough
- Death
- Diarrhea
- Dizziness
- Embolism
- Esophageal damage (including atrio-esophageal fistula)
- Fatigue
- Fever
- Headache
- Hemoptysis
- Hypotension/hypertension
- Infection (e.g. pericarditis, sepsis, urinary)
- Lightheadedness
- Myocardial infarction
- Nausea/vomiting
- Perforation
- Pericardial effusion
- Phrenic nerve injury
- Pleural effusion
- Pneumonia
- Pneumothorax
- Pseudoaneurysm
- Pulmonary edema
- Pulmonary hemorrhage
- Pulmonary vein dissection
- Pulmonary vein stenosis
- Shivering
- Shortness of breath
- Sore throat
- Transient ischemic attack
- Vagal nerve injury (e.g. gastroparesis)
- Vasovagal reaction
- Visual changes (e.g. blurred vision)

6 Instructions for use

6.1 Connecting the Arctic Front Advance Pro Cardiac Cryoablation Catheter

Use standard aseptic technique when removing the product from its sterile barrier and handling the catheter. During connection of the catheter, keep the connections dry at all times. Fluid incursion may lead to system malfunction. To connect the catheter, follow these steps. (For more detailed instructions, see the *CryoConsole Operator's Manual*.)

1. Connect the auto connection box to the CryoConsole.
Note: The ECG cable is not required for an Arctic Front Advance Pro Cryoballoon procedure, and should not be connected to the auto connection box.
2. Connect the catheter to a sterile electrical umbilical cable and then pass the other end of the cable out of the sterile field and connect it to the auto connection box.
3. Connect the catheter to a sterile coaxial umbilical cable and then pass the other end of the cable out of the sterile field and connect it to the CryoConsole.
4. Enable vacuum on CryoConsole.

6.2 Cryoablation

To use the catheter for a cryoablation procedure, follow these steps. (For more detailed instructions, see the *CryoConsole Operator's Manual*.)

Notes:

- Before introducing the catheter into the patient, test the deflection mechanism on the handle to ensure it is operational.
 - Always use the deflection mechanism on the handle to straighten the shaft before insertion or withdrawal of the catheter.
 - The Achieve family of mapping catheters is compatible for use with the Medtronic Arctic Front Advance Pro Cryoablation Catheter and may be used to support and position the catheter.
1. Using aseptic technique, create a vascular access with an appropriate introducer. Obtain left atrial transseptal access using a transseptal sheath, its dilator, and needle.
 - Place standard diagnostic pacing catheters.
 - Visualize left atrial anatomy to help select a balloon size. Select balloon size.
 - Selection of balloon size should be based on pulmonary vein diameter and shape, the surrounding anatomy, and desired position of the balloon outside the tubular portion of the pulmonary vein.
 2. Remove the transseptal sheath and dilator, leaving the guide wire positioned preferably in the left superior pulmonary vein.
 3. Advance the sheath and dilator over the wire into the left atrium.
 4. Slowly remove the guide wire and dilator from the sheath.
 5. Aspirate and flush the sheath.
 6. Obtain a catheter. Connect the Tuohy-Borst adapter and manifold to the push button Luer.
 7. Flush the guide wire lumen, ensuring there are no air bubbles in the guide wire lumen or connections.
 8. Load the circular mapping catheter or guide wire into the guide wire lumen.
 9. Pull the circular mapping catheter or guide wire back into the balloon catheter until it is just inside the balloon catheter tip.
 10. Flush the guide wire lumen again, ensuring there are no air bubbles in the guide wire lumen or connections.
 11. Retract the sleeve from the balloon onto the catheter shaft.
 12. Submerge and rinse the balloon to remove trapped air bubbles.
 13. Advance the sleeve back over the balloon while the balloon is still submerged.
 14. Place the sleeve adjacent to the hemostasis valve of the sheath.
Note: To avoid damaging the valve or introducing air, do not push the sleeve through the valve opening.
 15. Insert the catheter into the sheath.
Notes:
 - To avoid kinking the catheter, grip the catheter shaft at a distance no greater than 30 mm (1.2 in) from the sleeve during insertion.
 - Follow institutional procedures to aspirate and flush the sheath.
 16. Advance the balloon until the distal tip of the catheter aligns with the distal tip of the sheath using fluoroscopic guidance or other appropriate visualization techniques.
 17. Advance the circular mapping catheter or guide wire to the target pulmonary vein using fluoroscopic guidance or other appropriate visualization techniques.
 18. Advance the balloon over the circular mapping catheter or guide wire into the left atrium using fluoroscopic guidance or other appropriate visualization techniques.
 19. Inflate the balloon in the left atrium by pressing the Start button on the CryoConsole control panel for 2 s.
 20. Position the catheter at the ostium of the target pulmonary vein and not inside the tubular portion of the pulmonary vein.
 21. Verify the balloon position by injecting a mixture of 50/50 contrast/heparinized saline into the catheter guide wire lumen port or by using other appropriate techniques. Be sure to flush the guide wire lumen with heparinized saline after each contrast injection.

Notes:

- To improve balloon position and support, reposition the circular mapping catheter. If a stable balloon position cannot be obtained, exchange the circular mapping catheter for a guide wire.
- Before exchanging the circular mapping catheter for a guide wire, retract the catheter into the sheath.
- When using an auto injector for contrast delivery, ensure that the pressure limit does not exceed 500 psig.
- Follow contrast labeling and institutional procedures regarding the appropriate medical strategies to minimize the risk to the patient associated with using contrast.

22. Perform the cryoablation.**Notes:**

- During the initial ablation phase, monitor the balloon's position using appropriate visualization techniques. Moisture in the system may cause the inflated balloon to deflate and then re-inflate during the initial transition phases.
- Set the ablation duration on the CryoConsole screen.
- Physicians may modify the preset ablation duration based on clinical judgment. It is not recommended that the ablation duration exceeds 240 seconds.
- The balloon's outer diameter varies from inflation to cryoablation and may cause the balloon to shift.

23. Wait for the cryoablation phase to complete (at the end of the preset duration). The balloon remains inflated and the thawing phase begins.

Note: At any time, the ablation can be stopped by pressing the Stop Current Action button on the CryoConsole control panel or by pressing the Stop CryoAblation button on the screen.


24. During the thawing phase, observe the temperature indicator on the screen. The balloon deflates automatically when the temperature reaches 20°C.**25. Before moving the balloon, use appropriate techniques to ensure that the catheter is not adhered to tissue.****26. Determine effective ablation of the cardiac tissue by assessing electrical isolation of the pulmonary vein from the left atrium (entrance and exit block) after the cryoablation is complete.**

Note: As needed, perform additional treatments by positioning the balloon differently in the same pulmonary vein.

27. Position the catheter at the ostium of the next target pulmonary vein using the sheath, circular mapping catheter, or guide wire. Return to Step 17 and continue ablation.**28. To retract the balloon into the sheath, perform the following steps:**

Note: Ensure that the distal tip of the catheter is free to move to its maximum length.

- Use the deflection mechanism on the sheath handle to straighten the sheath. Use the deflection mechanism on the catheter handle to straighten the catheter.
- Inflate the balloon.
- Perform the following two steps simultaneously:
Advance the blue push button on the catheter handle, as shown in Figure 2. This causes the balloon to extend to maximum length and wrap tightly.
Deflate the balloon by pressing the Stop Current Action button on the CryoConsole control panel, or by selecting the Deflate Balloon option on the screen.
- Retract the catheter into the sheath.

Figure 2. Catheter handle**29. Remove the catheter from the patient.****30. After the procedure, follow the instructions to "Shut down the system" in the CryoConsole Operator's Manual. Follow the prompts to close the refrigerant tank and replace the cap on the coaxial port. Confirm these actions when prompted by the CryoConsole.****7 Specifications**

Catheter shaft outer diameter	3.5 mm (10.5 Fr; 0.14 in)
Tip length	8 mm (0.31 in)
Tip outer diameter	3.3 mm (9.9 Fr; 0.13 in)
Recommended introducer sheath	compatible Medtronic 12 Fr inner diameter sheath
Inner diameter of guide wire lumen	1.27 mm (0.05 in) nominal
Inflated balloon diameter	AFAPRO23 – 23 mm (0.91 in) AFAPRO28 – 28 mm (1.10 in)
Effective length (with balloon inflated)	95.0 ±2.0 cm (37.40 ±0.80 in)
Number of thermocouples	1
Environmental parameters	
Storage temperature	15°C to 30°C (59°F to 86°F)
Transit temperature	-35°C to 45°C (-31°F to 113°F); up to 85% relative humidity (non-condensing)
Operation	15°C to 30°C (59°F to 86°F) at altitudes less than 2400 m (8000 feet) above sea level

8 Medtronic limited warranty

For complete warranty information, see the accompanying limited warranty document.

9 Service

Medtronic employs highly trained representatives and engineers located throughout the world to serve you and, upon request, to provide training to qualified hospital personnel in the use of Medtronic products. Medtronic also maintains a professional staff to provide technical consultation to product users. For more information, contact your local Medtronic representative, or call or write Medtronic at the appropriate telephone number or address listed on the back cover.

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The screenshot shows the 'CryoTherapy' application interface. At the top, there's a header with the title 'CryoTherapy' and a small input field. Below it, a message says 'Touch a field to enter patient info. Touch "Therapy" to continue.' The main section is titled 'Patient Info' and contains several fields: Name, ID Number, Date of Birth, Gender, and Physician. The 'Date of Birth' field is currently selected, and a 'Date' selection window is overlaid on top. This window has a grid for selecting Year, Month, and Day. The Year is set to 1971, Month to 1, and Day to 11. There are 'OK' and 'Cancel' buttons at the bottom of the date selection window. At the very bottom of the screen, there is a navigation bar with buttons for 'Home', '1 Patient Info', '2 Therapy', and '3 Report'.

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Achieve Advance. Arctic Front. Medtronic

Explanation of symbols / Explicación des símbolos / Erläuterung der Symbole / Explicación de los símbolos / Verklaring van de symbolen / Spiegazione dei simboli / Symbolforklaring / Symboler / Επεξήγηση των συμβόλων / Symbolforklaring / Explicação dos símbolos / A szimbólumok magyarázata / Objaśnienie symboli / Vysvětlení symbolů / Объяснение символов / Vysvětlivky k symbolom / Sembollerin açıklaması / Explicação dos símbolos

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1. Description
The catheter is intended to be used with the Active Advance mapping catheter.
The mapping catheter is an intracardiac electrophysiology (IEP) recording catheter and can be used for cardiac stimulation during electrophysiology studies.
It is intended for use in the right atrium and right ventricle of the heart. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus.


Figure 1. Active Advance mapping catheter.
The mapping catheter should only be used with the Active Advance mapping catheter. It is not intended for use in the left atrium and left ventricle of the heart. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus.


Figure 2. Active Advance mapping catheter.
The mapping catheter should only be used with the Active Advance mapping catheter. It is not intended for use in the left atrium and left ventricle of the heart. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus. It is used to map the electrical activity of the heart and to determine the location of the arrhythmogenic focus.

1.1. Indications for use
The mapping catheter is used for the following purposes:
• Active Advance mapping catheter

1.2. Contraindications
The mapping catheter is contraindicated for the following purposes:
• Active Advance mapping catheter

1.3. Precautions for use
The mapping catheter should be used with the following precautions:
• Active Advance mapping catheter

1.4. Storage and handling
The mapping catheter should be stored and handled as follows:
• Active Advance mapping catheter

1.5. Maintenance
The mapping catheter should be maintained as follows:
• Active Advance mapping catheter

1.6. Disposal
The mapping catheter should be disposed of as follows:
• Active Advance mapping catheter

1.7. Packaging
The mapping catheter should be packaged as follows:
• Active Advance mapping catheter

1.8. Labeling
The mapping catheter should be labeled as follows:
• Active Advance mapping catheter

1.9. Instructions for use
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.10. Appendix
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.11. Bibliography
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.12. Glossary
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.13. Index
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.14. Annex
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.15. Declaration of conformity
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.16. Certificate of approval
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.17. Declaration of conformity
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.18. Certificate of approval
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.19. Declaration of conformity
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.20. Certificate of approval
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.21. Declaration of conformity
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.22. Certificate of approval
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.23. Declaration of conformity
The mapping catheter should be used as follows:
• Active Advance mapping catheter

1.24. Certificate of approval
The mapping catheter should be used as follows:
• Active Advance mapping catheter

Notes:

- Refer to the appropriate pacing and recording equipment operator manual for proper set up and operation.
- When connecting the catheter connecting cable (Model 2ACHC) to the EP recording system, pins 8 and 10 will not be used on the 8-electrode catheters (Models 2ACH15 and 2ACH20).
- Use of the Medtronic 8-pin cable (Model 930068) will result in a loss of functionality. It used with the 25 mm Achieve Advance mapping catheter (Model 2ACH25).

4. Advance the mapping catheter into the desired position.

Note: Mishandling of the mapping catheter could result in deformation of the distal loop. Warning: Always rotate the connector in a clockwise direction to position the mapping loop in the appropriate location to avoid tissue injury.

T Catheter specifications

Maximum catheter shaft size	3.2 Fr 1.0 mm (0.04 in)
Connector mating device	2-pin 1.0 mm (0.04 in)
Loop diameter	2ACH15: 15 mm (0.59 in) 2ACH20: 20 mm (0.79 in) 2ACH25: 25 mm (0.98 in)
Shaft length	Overall: 160 cm (52.8 in) Distal loop: 140 cm (45.9 in)
Number of electrodes	2ACH15: 8 2ACH20: 8 2ACH25: 10
Electrode length	1.0 mm (0.039 in)
Electrode Spacing	2ACH15: 4 mm (0.16 in) 2ACH20: 6 mm (0.24 in) 2ACH25: 8 mm (0.31 in)
Environmental parameters	
Operation	10 °C to 40 °C (50 °F to 104 °F)
Transit temperature	-20 °C to 50 °C (-4 °F to 122 °F), up to 500 minutes Humidity: 5% to 95%
Storage temperature	15 °C to 30 °C (59 °F to 86 °F)

8 Medtronic limited warranty

For catheter warranty information, see the accompanying limited warranty document.

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1 Description

Medvonic Confidential
Cathay Pacific

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As a result, the authors conclude that the use of a single, standardized, and validated instrument to assess the impact of the intervention on the target population is a necessary step in the development of a comprehensive evaluation of the intervention. The authors also note that the use of a single, standardized, and validated instrument to assess the impact of the intervention on the target population is a necessary step in the development of a comprehensive evaluation of the intervention.

[illegible]

the same time, the fact that the company is not a public company and is not subject to the same level of scrutiny as a public company is a significant advantage. The company's management team is experienced and has a strong track record of success. The company's financial performance is strong and consistent. The company's market position is strong and competitive. The company's growth potential is high. The company's risk profile is low. The company's valuation is reasonable. The company's management team is experienced and has a strong track record of success. The company's financial performance is strong and consistent. The company's market position is strong and competitive. The company's growth potential is high. The company's risk profile is low. The company's valuation is reasonable.

[illegible][illegible]

It is the policy of the Department of Education to ensure that all students have access to the same quality of education. This includes providing a safe and secure learning environment for all students. The Department is committed to ensuring that all students have access to the same quality of education, regardless of their background or circumstances. This includes providing a safe and secure learning environment for all students. The Department is committed to ensuring that all students have access to the same quality of education, regardless of their background or circumstances.

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5. Write the following:

© 1997 by The McGraw-Hill Companies, Inc. All rights reserved. Printed in the United States of America. This book is printed on acid-free paper.

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
- 1 AcNave Advanced Mapping-Karte

- 1 Achieve Advanced Shopping-Marketing
- Develop a Solid Business Plan

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1. The first step is to identify the problem. In this case, the problem is that the system is not working properly. The user has reported that the system is not working properly, and the user has provided some information about the problem. The first step is to identify the problem, and the user has provided some information about the problem. The first step is to identify the problem, and the user has provided some information about the problem.

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

Page 3, Line 1: Advance Advance handling page 31
 The above mentioned section is being requested in India as the above mentioned
 work is to be carried out by the above mentioned firm. The above mentioned firm
 is engaged in the above mentioned work. The above mentioned firm is engaged in the
 above mentioned work. The above mentioned firm is engaged in the above mentioned
 work. The above mentioned firm is engaged in the above mentioned work.

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CONCLUSIONS

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Medtronic has validated the representative electrode array and has not validated the individual electrodes. The representative electrode array is used for all testing and validation. The individual electrodes are used for specific testing and validation. The representative electrode array is used for all testing and validation. The individual electrodes are used for specific testing and validation.

Medtronic bagittal gait	
Forward temperature	15 °C (59 °F) to 35 °C (95 °F) (30 °F to 95 °F)
Transpiration rate	10 °C (50 °F) to 35 °C (95 °F) (30 °F to 95 °F)
Electrode array	24CH15, 6 mm (0.24 in) 24CH20, 8 mm (0.31 in) 24CH25, 10 mm (0.39 in) 24CH30, 12 mm (0.47 in) 24CH35, 15 mm (0.59 in) 24CH40, 18 mm (0.71 in) 24CH45, 20 mm (0.79 in) 24CH50, 25 mm (0.98 in) 24CH55, 30 mm (1.18 in) 24CH60, 35 mm (1.38 in) 24CH65, 40 mm (1.57 in) 24CH70, 45 mm (1.77 in) 24CH75, 50 mm (1.97 in) 24CH80, 55 mm (2.17 in) 24CH85, 60 mm (2.36 in) 24CH90, 65 mm (2.56 in) 24CH95, 70 mm (2.76 in) 24CH100, 75 mm (2.95 in) 24CH105, 80 mm (3.15 in) 24CH110, 85 mm (3.35 in) 24CH115, 90 mm (3.54 in) 24CH120, 95 mm (3.74 in) 24CH125, 100 mm (3.94 in) 24CH130, 105 mm (4.13 in) 24CH135, 110 mm (4.33 in) 24CH140, 115 mm (4.53 in) 24CH145, 120 mm (4.72 in) 24CH150, 125 mm (4.92 in) 24CH155, 130 mm (5.12 in) 24CH160, 135 mm (5.31 in) 24CH165, 140 mm (5.51 in) 24CH170, 145 mm (5.71 in) 24CH175, 150 mm (5.91 in) 24CH180, 155 mm (6.10 in) 24CH185, 160 mm (6.30 in) 24CH190, 165 mm (6.50 in) 24CH195, 170 mm (6.70 in) 24CH200, 175 mm (6.90 in) 24CH205, 180 mm (7.10 in) 24CH210, 185 mm (7.30 in) 24CH215, 190 mm (7.50 in) 24CH220, 195 mm (7.70 in) 24CH225, 200 mm (7.90 in) 24CH230, 205 mm (8.10 in) 24CH235, 210 mm (8.30 in) 24CH240, 215 mm (8.50 in) 24CH245, 220 mm (8.70 in) 24CH250, 225 mm (8.90 in) 24CH255, 230 mm (9.10 in) 24CH260, 235 mm (9.30 in) 24CH265, 240 mm (9.50 in) 24CH270, 245 mm (9.70 in) 24CH275, 250 mm (9.90 in) 24CH280, 255 mm (10.10 in) 24CH285, 260 mm (10.30 in) 24CH290, 265 mm (10.50 in) 24CH295, 270 mm (10.70 in) 24CH300, 275 mm (10.90 in) 24CH305, 280 mm (11.10 in) 24CH310, 285 mm (11.30 in) 24CH315, 290 mm (11.50 in) 24CH320, 295 mm (11.70 in) 24CH325, 300 mm (11.90 in) 24CH330, 305 mm (12.10 in) 24CH335, 310 mm (12.30 in) 24CH340, 315 mm (12.50 in) 24CH345, 320 mm (12.70 in) 24CH350, 325 mm (12.90 in) 24CH355, 330 mm (13.10 in) 24CH360, 335 mm (13.30 in) 24CH365, 340 mm (13.50 in) 24CH370, 345 mm (13.70 in) 24CH375, 350 mm (13.90 in) 24CH380, 355 mm (14.10 in) 24CH385, 360 mm (14.30 in) 24CH390, 365 mm (14.50 in) 24CH395, 370 mm (14.70 in) 24CH400, 375 mm (14.90 in) 24CH405, 380 mm (15.10 in) 24CH410, 385 mm (15.30 in) 24CH415, 390 mm (15.50 in) 24CH420, 395 mm (15.70 in) 24CH425, 400 mm (15.90 in) 24CH430, 405 mm (16.10 in) 24CH435, 410 mm (16.30 in) 24CH440, 415 mm (16.50 in) 24CH445, 420 mm (16.70 in) 24CH450, 425 mm (16.90 in) 24CH455, 430 mm (17.10 in) 24CH460, 435 mm (17.30 in) 24CH465, 440 mm (17.50 in) 24CH470, 445 mm (17.70 in) 24CH475, 450 mm (17.90 in) 24CH480, 455 mm (18.10 in) 24CH485, 460 mm (18.30 in) 24CH490, 465 mm (18.50 in) 24CH495, 470 mm (18.70 in) 24CH500, 475 mm (18.90 in) 24CH505, 480 mm (19.10 in) 24CH510, 485 mm (19.30 in) 24CH515, 490 mm (19.50 in) 24CH520, 495 mm (19.70 in) 24CH525, 500 mm (19.90 in) 24CH530, 505 mm (20.10 in) 24CH535, 510 mm (20.30 in) 24CH540, 515 mm (20.50 in) 24CH545, 520 mm (20.70 in) 24CH550, 525 mm (20.90 in) 24CH555, 530 mm (21.10 in) 24CH560, 535 mm (21.30 in) 24CH565, 540 mm (21.50 in) 24CH570, 545 mm (21.70 in) 24CH575, 550 mm (21.90 in) 24CH580, 555 mm (22.10 in) 24CH585, 560 mm (22.30 in) 24CH590, 565 mm (22.50 in) 24CH595, 570 mm (22.70 in) 24CH600, 575 mm (22.90 in) 24CH605, 580 mm (23.10 in) 24CH610, 585 mm (23.30 in) 24CH615, 590 mm (23.50 in) 24CH620, 595 mm (23.70 in) 24CH625, 600 mm (23.90 in) 24CH630, 605 mm (24.10 in) 24CH635, 610 mm (24.30 in) 24CH640, 615 mm (24.50 in) 24CH645, 620 mm (24.70 in) 24CH650, 625 mm (24.90 in) 24CH655, 630 mm (25.10 in) 24CH660, 635 mm (25.30 in) 24CH665, 640 mm (25.50 in) 24CH670, 645 mm (25.70 in) 24CH675, 650 mm (25.90 in) 24CH680, 655 mm (26.10 in) 24CH685, 660 mm (26.30 in) 24CH690, 665 mm (26.50 in) 24CH695, 670 mm (26.70 in) 24CH700, 675 mm (26.90 in) 24CH705, 680 mm (27.10 in) 24CH710, 685 mm (27.30 in) 24CH715, 690 mm (27.50 in) 24CH720, 695 mm (27.70 in) 24CH725, 700 mm (27.90 in) 24CH730, 705 mm (28.10 in) 24CH735, 710 mm (28.30 in) 24CH740, 715 mm (28.50 in) 24CH745, 720 mm (28.70 in) 24CH750, 725 mm (28.90 in) 24CH755, 730 mm (29.10 in) 24CH760, 735 mm (29.30 in) 24CH765, 740 mm (29.50 in) 24CH770, 745 mm (29.70 in) 24CH775, 750 mm (29.90 in) 24CH780, 755 mm (30.10 in) 24CH785, 760 mm (30.30 in) 24CH790, 765 mm (30.50 in) 24CH795, 770 mm (30.70 in) 24CH800, 775 mm (30.90 in) 24CH805, 780 mm (31.10 in) 24CH810, 785 mm (31.30 in) 24CH815, 790 mm (31.50 in) 24CH820, 795 mm (31.70 in) 24CH825, 800 mm (31.90 in) 24CH830, 805 mm (32.10 in) 24CH835, 810 mm (32.30 in) 24CH840, 815 mm (32.50 in) 24CH845, 820 mm (32.70 in) 24CH850, 825 mm (32.90 in) 24CH855, 830 mm (33.10 in) 24CH860, 835 mm (33.30 in) 24CH865, 840 mm (33.50 in) 24CH870, 845 mm (33.70 in) 24CH875, 850 mm (33.90 in) 24CH880, 855 mm (34.10 in) 24CH885, 860 mm (34.30 in) 24CH890, 865 mm (34.50 in) 24CH895, 870 mm (34.70 in) 24CH900, 875 mm (34.90 in) 24CH905, 880 mm (35.10 in) 24CH910, 885 mm (35.30 in) 24CH915, 890 mm (35.50 in) 24CH920, 895 mm (35.70 in) 24CH925, 900 mm (35.90 in) 24CH930, 905 mm (36.10 in) 24CH935, 910 mm (36.30 in) 24CH940, 915 mm (36.50 in) 24CH945, 920 mm (36.70 in) 24CH950, 925 mm (36.90 in) 24CH955, 930 mm (37.10 in) 24CH960, 935 mm (37.30 in) 24CH965, 940 mm (37.50 in) 24CH970, 945 mm (37.70 in) 24CH975, 950 mm (37.90 in) 24CH980, 955 mm (38.10 in) 24CH985, 960 mm (38.30 in) 24CH990, 965 mm (38.50 in) 24CH995, 970 mm (38.70 in) 24CH1000, 975 mm (38.90 in)

2. For the electrode array, the representative electrode array is used for all testing and validation. The individual electrodes are used for specific testing and validation. The representative electrode array is used for all testing and validation. The individual electrodes are used for specific testing and validation.

M4265008001 Rev 1A

[illegible]

From 1982 to 1990, the number of people in the United States who were employed in the health care industry increased by 150 percent. This increase was due to a number of factors, including the aging of the population, the increasing prevalence of chronic diseases, and the growing demand for health care services. The health care industry is now one of the largest and fastest-growing sectors of the U.S. economy.

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The following table shows the results of the regression analysis for the dependent variable $\ln(\text{Sales})$. The independent variables are $\ln(\text{Price})$, $\ln(\text{Advertising})$, and $\ln(\text{Market Size})$. The table includes the coefficient estimates, standard errors, t-statistics, and p-values for each variable.

12 **strongly supports** **strongly opposes** **neutral**

4) **Strongly supports Federal Authority** (strongly by

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MAZ

Toplak shladovani
15 °C ± 30 °C (59 °F ± 88 °F)

Chemical name	Chemical structure	Chemical formula	Molecular weight	Boiling point (°C)	Melting point (°C)	Density (g/cm³)	Refractive index	Log P	Water solubility (mg/L)	Biodegradability	Environmental fate	Ecotoxicity	Human toxicity	Regulatory status	Other information
1,1,1-Trichloroethane		<chem>CCl3</chem>	131.4	74	-136	1.29	1.29	2.23	130	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,2-Trichloroethane		<chem>CCl2CCl</chem>	131.4	112	-134	1.29	1.29	2.23	130	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1-Dichloroethane		<chem>CCl2</chem>	98.96	36	-98	1.20	1.20	1.98	100	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2-Dichloroethane		<chem>CClCCl</chem>	98.96	83	-100	1.25	1.25	1.98	100	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,1,2-Tetrachloroethane		<chem>CCl3CCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,2,2-Tetrachloroethane		<chem>CCl2CCl2</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,2,3-Tetrachloroethane		<chem>CCl2CClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,3,3-Tetrachloroethane		<chem>CCl2CClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,3,4-Tetrachloroethane		<chem>CCl2CClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,1,4,4-Tetrachloroethane		<chem>CCl2CClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,3,4-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,3,5-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,4,4-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,4,5-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,5,5-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,3,3,4-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,3,4,4-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,3,4,5-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,3,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,2,6,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,3,6,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,4,4,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,4,5,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,5,5,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent
1,6,6,6-Tetrachloroethane		<chem>CClCClCClCCl</chem>	167.8	106	-136	1.47	1.47	2.68	10	Yes	Atmospheric degradation	Low	Low	Not regulated	Used as a solvent

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

1.

Abstract The purpose of this study was to determine the effect of a 10-week, 1000 kcal energy deficit diet on the body composition and physical fitness of obese women. The subjects were 20 obese women who were randomly assigned to either a diet or a control group. The diet group was instructed to consume a diet that was 1000 kcal less than their maintenance level, while the control group was instructed to consume a diet that was equal to their maintenance level. The subjects were monitored for 10 weeks, and their body composition and physical fitness were measured at the beginning and end of the study. The diet group showed a significant decrease in body weight, body fat percentage, and waist circumference, while the control group showed no significant changes. The diet group also showed a significant increase in physical fitness, while the control group showed no significant changes. The results of this study suggest that a 10-week, 1000 kcal energy deficit diet is effective in improving body composition and physical fitness in obese women.

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© 1999 by the American Psychological Association. All rights reserved. This article is intended solely for the personal use of the individual user and is not to be disseminated broadly.

[illegible]

Abstract 2. Study's main purpose is to help Achieve Agenda

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4 ~~tratat~~ **tratat** ~~an~~ **an** ~~si~~ **si** ~~pre~~ **pre** ~~ventiv~~ **ventiv** ~~ne~~ **ne** ~~op~~ **op** ~~ar~~ **ar** ~~in~~ **in** ~~la~~ **la**

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Medtronic Confidential
CathAF_iFU_R1.5

Medtronic Confidential
CathAF_iFU_R1.5

Medtronic Confidential

Refer to the "CathAF_iFU" category
in doc#162256 for Printing
Instructions.

M4265000, CathAF, Rev 1.0, 1.00 mm
4" x 12" (100 mm x 300 mm)

Medtronic Confidential
CathAF_JFU_R1.5

M4265000001 Rev 1.0

Refer to the "CathAF_JFU" category
in doc#163256 for Printing
Instructions.

MA26500B001 Rev 1A
4" x 10" x 100 mm x 200 mm

Medtronic Confidential
CathAF_IFU_R1.5

Medtronic

Manufacture:
710 Medtronic Parkway
Minneapolis, MN 55432
USA
www.medtronic.com
+1 763 514 4000

ex hcr

Authorized Representative in the European Community

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Australia
Medtronic Australasia Pty Ltd
97 Walerloo Road
North Ryde, NSW 2113
Australia

Canada
Medtronic of Canada Ltd
99 Hereford Street
Brampton, Ontario L6Y 0R3
Canada
+1 905 460 3800

Technical manual
www.medtronic.com/medtronic

© Medtronic, Inc. 2010
MA26500B001 1A
2010-03-16

* M 4 2 6 5 0 0 B 0 0 1 *

MA26500B001 Rev 1A

Refer to the "CathAF_IFU" category
in doc#163258 for Pinning
Instructions

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

Declaration Letter


Medtronic Inc., a corporation organized under the laws of the United States, having its principal place of business at 710 Medtronic Parkway, Minneapolis, MN 55432, USA, hereby swear (or affirm) that the attached documentation is true and unaltered copies of the original document.

Donna J. Bergin
Principal Regulatory Affairs Specialist

3-28-2018
Date

Notary stamp

3/28/2018

As a registered Notary Public for the State of Minnesota,
I hereby certify that this is a true and accurate copy of the
original document. The original is filed at Medtronic, Inc.

Medtronic

2ACHC

Catheter connecting cable
Câble de connexion de cathéter
Katheteranschlusskabel
Cable de conexión de catéter
Katheter-verbindingskabel
Cavo di collegamento per cateteri
Forbindelseskabel til kateter
Kateteranslutningskabel
Καλώδιο σύνδεσης καθετήρα
Kateterforbindelseskabel
Cabo de ligação de cateteres
Katéter-csatlakoztatókábel
Przewód łączący cewnik
Spojovací kabel katétru
Соединительный кабель катетера
Pripájací kábel katétra
Kateter bağlantı kablosu
Cabo de conexão de cateter

Phuong

3/28/2018

Technical Manual • Manuel technique • Technisches Handbuch • Manual técnico • Technische handleiding • Manuale tecnico • Teknisk håndbok • Teknisk manual • Τεχνικό εγχειρίδιο • Teknisk håndbog • Manual técnico • Műszaki leírás • Instrukcja techniczna • Technická příručka • Техническое руководство • Technická příručka • Teknik El Kitabı • Manual técnico

! USA

R_x only

CE0088

2015

The following list includes trademarks or registered trademarks of Medtronic in the United States and possibly in other countries. All other trademarks are the property of their respective owners.

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1 Description

The sterile, single use only Medtronic catheter connecting cable, Model 2ACHC, provides conduction from the proximal end of the mapping catheter to standard shielded ECG pins that connect into standard EP recording and pacing equipment. The cable is not intended for reuse or resterilization.

2 Intended use

The cable is designed for use only with the Achieve family of mapping catheters. It is intended to provide connection of the catheter to a standard ECG interface box.

3 Warnings and precautions

For single use only – This cable is intended only to be used once for a single patient. Do not reuse, reprocess, or resterilize these devices for purpose of reuse. Reuse, reprocessing, or resterilization may compromise the structural integrity of the device or create a risk of contamination of the device that could result in patient injury, illness, or death.

Do not resterilize – Do not resterilize the cable for purpose of reuse. Resterilization may compromise the structural integrity of the device or create a risk of contamination from the device that could result in patient injury, illness, or death.

Device disposal – Dispose of the cable according to hospital biohazard requirements.

Expiration date – Check to verify the cable is within the expiration date. Do not use if the product date has expired.

Sterile package inspection – Inspect the sterile packaging prior to use. If the sterile packaging of the device is damaged, do not use. Contact your local Medtronic representative.

Cable integrity – Do not use the cable if it is kinked or damaged. If the cable becomes kinked or damaged while in use, remove it and use a new cable.

4 Instructions for use

1. Connect the shielded tip jacks from the cable to the ECG interface box.

Note: When connecting to the EP recording system, electrodes 9 and 10 will not be used for models within the Achieve family that have 8 electrodes.

2. Connect the cable single connector to the mapping catheter electrical connector.

4.1 Cable specifications

Length	1.96 m (77.0 in)
Environmental parameters	
Operation	10 °C to 40 °C (50 °F to 104 °F)
Transit temperature	-40 °C to 70 °C (-40 °F to 158 °F); 95% relative humidity (non-condensing)
Storage temperature	15 °C to 30 °C (59 °F to 86 °F)

5 Medtronic disclaimer of warranty

For cable disclaimer of warranty, see the disclaimer of warranty packaged with the cable.

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Circular mapping catheter compatibility – Use only Medtronic circular mapping catheters compatible with the inner lumen of the Arctic Front Advance Pro Cryoballoon.

медицинскими
изделиями / The
possibility and
ways of
integrating with
other medical
device

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Guide wire compatibility – Use only 0.081 cm (0.032 in) or 0.089 cm (0.035 in) guide wires with the catheter.

Steerable sheath compatibility – Use only a compatible Medtronic 12 Fr inner diameter sheath with the Arctic Front Advance Pro Cryoballoon.

System compatibility – Use only Medtronic cryoablation catheters, refrigerant tanks, and components with the CryoConsole.

The catheter shall directly interface with the following market released Medtronic devices and/or accessories as per Table:

Table: Catheter interfacing with market released Medtronic devices and/or accessories.

CFN	Description
4FC12	12Fr FlexCath Advance™ Steerable Sheath
990063-015	Achieve™ Mapping Catheter (15mm)
990063-020	Achieve™ Mapping Catheter (20mm)
2ACH15	Achieve™ Advance Mapping Catheter (15mm)
2ACH20	Achieve™ Advance Mapping Catheter (20mm)
2ACH25	Achieve™ Advance Mapping Catheter (25mm)
2035U, 2035UC	Sterile Electrical Umbilical
203CX, 203CXC	Sterile Coaxial Umbilical
20MRK	Manual Retraction Kit
AFAKITU	Cryoballoon Procedure Kit

The catheter shall interact with, but is not directly connected to, the following market released Medtronic devices and/or accessories as per Table:

Table: Catheter interaction with other market released Medtronic devices and/or accessories.

CFN	Description
209F3	Freezor ® MAX 3
209F5	Freezor ® MAX 3
2037A	Auto Connection Box
106E2	CrvoConsole Gen V Universal (model # 10000-008-02)
106A2-K	CrvoConsole Gen V Universal (model # 10000-008-03)
106A3	CrvoConsole Gen V Universal (model # 10000-008-04)
106E2 China	CrvoConsole Gen V Universal (model # 10000-008-05)
106E2	CrvoConsole Gen V RoHS (model # 10000-008-06)
106A2-K	CrvoConsole Gen V RoHS (model # 10000-008-07)
106A3	CrvoConsole Gen V RoHS (model # 10000-008-08)
1038 series	Power Cords
1035 series	Scavenging Hoses
1036 series	Adapters

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




Connecting Cable Model 2ACHC
Achieve Catheter Connecting Cable, Model 2ACHC

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


Item	Description
1	Cable, grey
2	Connector, redel, plug, in-line, 14 pins, black
3	Conn, 2 mm
4	Mold label Insert
5	Kester lead free solder
6	Lead wire, with PVC insulation
7	Wrap around label (UDI label)
8	Hot glue for potting
9	Black shroud material
10	Inner yolk material
11	Outer yolk material
12	Outer yolk colorant (black)
13	Kester flux
14	Heat shrink tubing
15	Splice band
16	Adhesive for rebel threads -- (Clear)

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

2ACHC

Catheter Connecting Cable

Package contents: Cable Product documentation

CS: Spojovací kabel katétru, DA: Kateterforbindelseskabel, DE: Katheteranschlusskabel,
EL: Καλώδιο σύνδεσης καθετήρα, ES: Cable de conexión de catéter, FR: Câble de connexion de cathéter,
HU: Katéter-csatlakoztatókábel, IT: Cavo di collegamento per cateteri, NL: Katheter-verbindingskabel,
NO: Forbindelseskabel til kateter, PL: Przewód łączący cewnik, PT: Cabo de conexão de cateter,
P2: Cabo de ligação de cateteres, RU: Соединительный кабель катетера, SK: Pripájací kábel katétra,
SV: Kateteranslutningskabel, TR: Kateter Bağlantı Kablosu

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medical products and their features

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the best course of treatment for a patient. Electrophysiology catheters typically consist of a polymer shaft with multiple recording electrodes mounted on the distal end. These electrodes are connected to fine wires that run inside of the shaft and through an attached manipulating handle. The wires terminate in a connector that can be connected to electrophysiology recording and cardiac pacing equipment. The Achieve Advance Mapping Catheter includes these technical characteristics as well.

The catheters are inserted into the patient's vasculature, typically through an introducer sheath via the femoral vein. The catheters are advanced into the right atrium. If access to the left atrium is desired, a transseptal sheath is used or deployed by advancing the sheath over the catheter and through the septum. The catheter is then advanced into the left atrium through the sheath.

The Achieve Advance Mapping Catheter is navigated to the target cardiac site. Once at the target site, the Achieve Advance Mapping Catheter maps using the multiple electrodes surrounding the distal ring or loop, oriented in a plane at approximate right angles to the main body of the catheter. When used with the Arctic Front family of cryoballoon catheters, the cryoballoon catheter is passed through the FlexCath sheath lumen over an Achieve Advance Mapping Catheter and into the left atrium in an uninflated state. Refer to Figure. The use of the Achieve Advance Mapping Catheter offers the following:

- Eliminates the need for a second transseptal puncture
- Minimizes the number of catheter exchanges through the FlexCath® sheath
- Measures PV potentials pre-, intra-, and post-cryoablation



Figure: Use of the Achieve Advance Mapping Catheter in Cryoballoon Procedure

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



Vladimir
Regulatory Affairs Specialist

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

Vienna Greg
Regulatory Affairs Specialist

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


Vladimir Gregor
Regulatory Affairs Specialist

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




ACHIEVE ADVANCE™ 2ACH20
MAPPING CATHETER

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




Vladimir Gregorovich
Regulatory Affairs Specialist

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ACHIEVE ADVANCE™ 2ACH25 MAPPING CATHETER

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

Characteristics	Mapping Catheter ACHIEVE ADVANCE.	Mapping Catheter ACHIEVE ADVANCE.	Mapping Catheter ACHIEVE ADVANCE.
-----------------	--	--	--

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Date of production / Series / Expiration date: see packaging	See Instructions See symbols on the package
Not toxic Pyrogen-free	

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Characteristics	
Characteristics	Meaning
Maximum catheter shaft size	3.3 Fr 1.09 mm (0.043 in)
Compatible mating device minimum internal diameter	3.4 Fr 1.12 mm (0.044 in)
Loop diameter	2ACH15: 15 mm (0.59 in) (+2.0/-1.0 mm) 2ACH20: 20 mm (0.79 in) (+2.0/-1.0 mm) 2ACH25: 25 mm (0.98 in) (+2.0/-1.0 mm)
Shaft length	Overall: 165 cm (65.0 in) Effective: 146 cm (57.5 in)
Number of electrodes	2ACH15: 8 2ACH20: 8 2ACH25: 10
Electrode length	1,0414 mm (0,041 in) ± 0,0762 mm (0,003 in)
Length of the introducer	139,7 ± 6,35 mm
Electrode Spacing	2ACH15: 4,1 ± 0,5 mm 2ACH20: 5,8 ± 0,5 mm 2ACH25: 5,8 ± 0,5 mm
Minimum tensile strength, N	Electrode to distal tubing: 5N (1,13ibf) Distal tubing to proximal joint: 22N (5,0 ibf) Electrical connector shell to proximal shaft: 44N (1,13ibf) Ball tip to distal tubing: 5N (1,13ibf)

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
Notes:

Two people are required for this procedure.

Unless instructed otherwise by a system notice, do not disconnect the coaxial and electrical umbilical cables until indicated in these steps.


Prepare the equipment

1. Unpack and assemble the manual retraction kit.
2. Verify that the stopcock and syringe are assembled as shown.


Release the vacuum and attach the coaxial connector

1. On the Therapy tab, touch **Release Vacuum** on the screen.
2. Disconnect the electrical umbilical cable and the coaxial umbilical cable from the catheter, as shown:

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


3. Push the syringe plunger in.
4. Connect the manual retraction kit coaxial connector to the catheter coaxial connector:


Apply vacuum to aspirate the catheter


1. Rotate the stopcock handle to close flow to the vent:

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

2. Pull the plunger out of the syringe to approximately 60 cc.
3. Rotate the plunger to lock the syringe:


4. Rotate the stopcock handle to close flow to the catheter inlet:

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


5. Disconnect the syringe from the stopcock.
6. Rotate the plunger to unlock it.
7. Push in the syringe plunger to empty the syringe.
8. Reconnect the syringe to the stopcock.
9. Apply a second syringe stroke of vacuum to the catheter by repeating steps 1 through 8.

Prepare and perform saline inflation of the balloon


1. Ensure the flow to the catheter inlet is closed:


2. Disconnect the syringe from the stopcock:

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


3. Fill the syringe with 30-50 cc of saline.
Note: Minimum quantity is 20 cc.


4. Connect the syringe to the stopcock:

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

5. Rotate the stopcock handle to close flow to the vent:


6. Inject the saline into the catheter (minimum 20 cc):

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


Tension and deflate the balloon

1. Retract the saline, using two people, as follows:


- Advance the blue push button on the catheter handle.
- At the same time, pull on the plunger to retract the saline.

2. Rotate the stopcock handle to close flow to the catheter inlet:

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3. Consult the technical manual for the catheter for steps to remove the catheter from the patient.

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The manual retraction kit provides an alternative means of inflating and deflating the balloon catheter, in the unlikely case that the normal methods of inflating and deflating the balloon cannot be used. The manual retraction kit is a sterile part intended for single use.




Figure A-1. Manual Retraction Kit

1. Coaxial connector
2. Stopcock
3. Syringe
4. Stopcock handle
5. Closed position indicator
6. Syringe connector
7. Coaxial connector inlet
8. Vent
9. Locking pin



The manual retraction kit has three main parts:

- **Coaxial connector** — connects the manual retraction kit to the coaxial connector of the catheter handle.
- **Stopcock** — a three-position valve that allows and restricts flow. The closed position indicator on the stopcock handle points to the closed passage for each stopcock position. (The stopcock is shipped preattached to the coaxial connector.)
- **Syringe** — creates vacuum to aspirate the catheter; injects saline solution into the balloon during inflation; extracts saline to collapse balloon.

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
Item	Description
1	Coaxial vacuum tube
2	Male coaxial fitting RoHS
3	Male luer to barb fitting
4	Snap luer lock ring
5	Vacuum o-ring
6	Injection o-ring
7	Male coaxial connector RoHS
8	60 CC vaclock syringe
9	THF solvent
10	Loctite 4307
11	High pressure stopcock

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Manual retraction kit is to be used with CryoBalloon catheter only to deflate them. The sterile 20MRK Manual Retraction Kit is designed for use with compatible Medtronic cryoablation balloon catheters.

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
Manual Retraction Kit, model 20MRK:

Contents of delivery:
- coaxial connector;
- stopcock;
- syringe.
- Product documentation.

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

Vladia Gregorski
Regulatory Affairs Specialist

Medtronic
 Medtronic CryoCath LP, Canada
 9000 Autoroute Transcanadienne,
 Pointe-Claire, Quebec H9R 5Z8, Canada
 www.medtronic.com

Vladia Gregorski
 Regulatory Affairs Specialist

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Characteristics primary packaging:

Characteristics	Manual Retraction Kit, model 20MRK
Primary packaging mass, g	72,4
Packaging dimensions (Length x Width), mm	302,5 x 204

* The tolerance is $\pm 10\%$, unless otherwise specified.

Characteristics secondary packaging:

Characteristics	Manual Retraction Kit, model 20MRK
Secondary packaging mass, g	210,6
Packaging dimensions (Length x Width x Height), mm	320 x 210 x 38,5

* The tolerance is $\pm 10\%$, unless otherwise specified.

Information on the Russian-language sticker
For sale on the territory of the Russian Federation, the final packaging is accompanied by an additional sticker.

Products come with a sticker that can contain the following information:

- The name of the medical device in Russian;
- Name of the manufacturer's factory;
- Country of origin;
- UPN code (unique product code);
- CFN code (model number);
- Date and number of the Registration Certificate;
- Mark of conformity in Russia;
- Warning and informative signs and inscriptions (if applicable).

Sample sticker:

	Name of medical device	
	Manufacturer:	Code UPN
	Country of origin:	Code CFN
	Manufacturer representative in Russia:	
	Registration Certificate № of	
	Date of production / Series /	See Instructions
	Expiration date: see packaging	See symbols on the package
	Not toxic Pyrogen-free	

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