
Medtronic

Declaration of Authenticity


Medtronic Sofamor Danek USA, Inc. located at 1800 Pyramid Place Memphis, TN 38132 USA, hereby declares that the IFU in the Appendix to Operating document, attached herein, to be a true and accurate copy of the original kept on file at Medtronic Sofamor Danek USA, Inc.


Signature: 
Rong Guo
Principle Regulatory Affairs Specialist

Date: 22 Feb 2023

COUNTY OF SHELBY
STATE OF TENNESSEE

SUBSCRIBED AND SWORN TO ME BEFORE THIS 22TH DAY OF FEBRUARY 2023.

Name

FOR INTERNATIONAL AUDIENCES ONLY

US audiences refer to www.medtronic.com/manuals

Medtronic Reusable Instruments and Accessories


US Audiences STOP - refer to www.medtronic.com/manuals

0380035 Rev. F

2020-04-08

Medtronic

* 0 3 8 0 0 3 5 - F *

Electronic instructions for use/manuals

To view, download, print, or order manuals for this product, go to www.medtronic.com/manuals, or contact a Medtronic representative.

Software requirements

Manuals can be viewed using a current version of any major internet browser. For best results, use Adobe Acrobat® Reader with the browser.

Paper manuals

Paper manuals are available to customers free of charge. They should arrive in 3-7 days. To order, go to www.medtronic.com/manuals, or contact a Medtronic representative.

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Physicians should review manuals before implanting, using a device, or performing a follow-up session.

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IMPORTANT INFORMATION FOR MEDTRONIC REUSABLE INSTRUMENTS AND ACCESSORIES

DESCRIPTION

Reusable instruments are manufactured from a variety of materials commonly used in orthopedic and neurological procedures, which meet applicable national and/or international standards. These instruments may be shipped in accessories.

Accessories include cases, trays, lids, caddies, and modules fabricated from a variety of materials, which meet national and international standards, and are commonly used to enclose, protect, and organize Medtronic orthopedic and neurological reusable instruments.

INTENDED USE

Reusable instruments are intended for use in orthopedic and neurological surgical procedures to manipulate tissue, bone, or for use with other devices in surgery. An instrument may incorporate a measuring function, which has uses as described on the label and the instrument.

Note: do not implant instruments.

Accessories are intended for use in healthcare facilities to organize, enclose, sterilize, transport, and store medical devices and other instrumentation between surgical and other medical uses. Accessories are not intended on their own to maintain sterility; they are intended to be used in conjunction with a legally marketed and/or appropriate sterilization wrap.

Medtronic does not warrant use of instruments, accessories, or any component parts upon which repairs have been made or attempted, except as performed by Medtronic or an authorized Medtronic repair representative. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

WARNINGS

- Breakage, slippage, misuse, or mishandling of instruments or accessories, such as on sharp edges, may cause injury to the patient, surgical personnel, or processing personnel.
- Improper maintenance, handling, or inadequate cleaning procedures can render instruments or accessories unsuitable for the intended purpose and may present a danger to the patient, surgical personnel, or processing personnel.
- There are additional risks associated with instruments used for bending and cutting rods. Use of these types of instruments can cause injury to the patient by virtue of the extremely high forces required. Do not cut rods in situ. Breakage of an instrument or implant in this situation could be extremely hazardous. Physical characteristics required for many instruments do not permit them to be manufactured from implantable materials. If broken fragments of instruments remain in the body of a patient, there could be allergic or infectious consequences.
- The surgeon should use extreme caution when working in close proximity to vital organs, nerves, or vessels. In addition, excessive force should not be used when positioning instruments or implants, as it could cause injury to the patient.
- Correct handling of accessories is extremely important. Do not modify accessories. Do not notch or bend accessories. Notches, scratches, wear, or other damage in accessories occurring during surgery may contribute to breakage.
- Immediate-Use Steam Sterilization is only intended for individual instruments and should only be considered under urgent clinical situations and when approved by local policies.

PRECAUTIONS

- Excessive force applied by instruments can dislodge implanted devices.
- Extreme care should be taken to ensure instruments and accessories remain in good working order. During the procedure, the proper functioning of instruments and accessories is extremely important. Instruments and accessories should not be bent or damaged. Misuse of instruments and accessories, corrosion, binding, scratching, loosening, bending, or fracture of components of instruments or accessories may inhibit or prevent proper function.
- Do not use instruments or accessories for any action for which they are not intended.
- Regularly review the operational state of all instruments and, if necessary, contact Medtronic or an authorized Medtronic repair representative. To avoid injury, the instrument should be carefully examined for functionality or damage prior to use. A damaged instrument should not be used. Additional back-up instruments should be available.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful use of instruments. Proper selection and compliance of the patient greatly affect results. Read and follow all product information supplied by the manufacturer of devices.
- Extra care should be taken when using instruments in pediatric patients since these patients can be more susceptible to stresses required in use of instruments.
- Ensure instruments or accessories with a measuring function are not worn and any surface markings are clearly visible.
- Personnel should wear all appropriate Personal Protective Equipment (PPE) as recommended by the sterilizer manufacturer.

When explanting and/or disposing of a device, be sure to avoid exposure to bodily substances such as blood, tissue, etc., as contact could lead to infection or disease. Always wear and use proper equipment, taking special care with sharp objects and needles. Follow your healthcare center's policy regarding both the disposal of devices and any events of exposure.

POSSIBLE ADVERSE EFFECTS

- Nerve damage, paralysis, pain, or damage to soft tissue, visceral organs, or joints.
- Infection if instruments or accessories are not properly cleaned and sterilized.
- Dural leak in cases of excessive load application.
- Impingement of close vessels, nerves, and organs by slippage or misplacement of the instrument.
- Damage due to spontaneous release of clamping devices or spring mechanisms of certain instruments.
- Cutting of skin or gloves of surgical or processing personnel.
- Bone fracture.
- Tissue damage to the patient, physical injury to surgical personnel, and/or increased operating time that may result from accidental disassembly of multi-component instruments occurring during surgery.
- Methods of use are determined by the user's experience and training in surgical procedures. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise results.

Physician note: although the physician is the learned intermediary between the company and patient, the important medical information in this document should be conveyed to the patient.

These devices should be used only by physicians familiar with the devices, their intended use, and available surgical techniques. In addition, accessories should be used only by operators familiar with the device(s), its intended use, and processing information.

PACKAGING

Instruments may be shipped in stand-alone packaging or in accessories. Stand-alone packaging, accessories, and contained medical devices should be intact upon receipt. Accessories and contained medical devices should be checked for completeness and signs of damage prior to use. Damaged packages and accessories should not be used and should be returned to Medtronic.

EXAMINATION

Instruments must always be examined pre/post cleaning and prior to surgery. Examination should be thorough and must include a visual and functional inspection of working surfaces, pivots, racks, spring or torsional operation, cleanliness of location holes or cannulations, and the presence of any cracks, bending, deformation, or distortion, and that all components are present. Never use instruments or accessories with obvious signs of excessive wear, damage, or that are incomplete or otherwise non-functional.

Accessories must be examined prior to use. Examination should be thorough and must include a visual and functional inspection. Examination must also ensure accessories are complete. Do not process damaged accessories.

Visual Inspection

Ensure the following:

- Laser markings, engravings, and other markings are legible.
- No cracks, deformation, or distortion are present.
- No cracks are present in instrument handles or any part of the instrument.
- No discoloration, corrosion, stains, or rust are present. If present, wipe the device clean in accordance with the instructions in the PROCESSING section.
- There is no handle/shaft separation, and the handle-to-shaft connection is secure.
- No cuts or gouges in the silicone of the brackets or instrument handles are present.
- No cracks or peeling of nylon coating on metal brackets.
- There is no damage (cuts, tears, etc.) to insulation.
- There is no damage to working ends or tips. The working end should be free of cracks, gouges, and other damage. When applicable, the working end should be sharp.
- There is no damage to instrument threads.
- All parts are present and free of damage and deterioration. Examples of parts that may be missing, loose, or damaged include set screws, springs, curved springs, pins, and prongs.
- Mating ends are free of damage (nicks, gouges, bends, etc.) that would interfere with the mating function.
- Cannulated instruments with a guide wire or other insertion tool are visually checked.

Functional Inspection

As applicable, ensure the following:

- Moving parts move freely, without sticking, binding, or grinding.
- Springs return the handle of the instrument to its original position.
- Retention tabs hold appropriate mating parts and are not damaged.
- The instrument will function as intended with appropriate mating parts.
- Ball detents will hold mating parts and are free from damage.
- Sharp edges are sharp to the touch and are not dull, have no nicks, or any other damage.
- Tips meet when appropriate.
- Ratcheting mechanisms are functional. This includes handles, latches, and other mechanisms. All teeth should be present and functional.
- Driver tips are not worn beyond functional use. If necessary, mate the instrument with the appropriate part.
- Latches, caddies, modules, and lids fit correctly.

If the device passes the visual and functional inspection as described in the IFU, the instrument can be expected to perform as intended for the perioperative period.

PROCESSING - GENERAL CONSIDERATIONS

Cleaning is the removal of organic soil. Effective cleaning:

- Avoids organic soil transfer from one patient to another.
- Allows for successful follow-up sterilization. Adequate processing is contingent upon thoroughness of cleaning.

Cleaning is the initial step. Sterilization occurs later in processing and is intended to kill microorganisms to reduce likelihood of transmission and possibility of infection. To ensure acceptable processing, there should be no delay between steps in this document.

Instruments may be safely and effectively processed using manual or a combination of manual and automated washer-disinfector (enzymatic or alkaline) cleaning procedures outlined in this document. Accessories may be safely and effectively processed using manual or a combination of manual and automate washer-disinfector (enzymatic) cleaning procedures outlined in this document.

The processor should comply with local laws and ordinances in countries where processing requirements are more stringent than those detailed in this document.

Bloodborne Pathogens

Universal precautions for handling these instruments after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030 regarding occupational exposure to bloodborne pathogens.

Cleaning Agents and Cleaning Tools

Instruments: use of neutral detergents, enzymatic cleaning agents (7.0–8.0pH), alkaline cleaning agents (8.0–11.0pH), soft bristled brushes, and soft pipe cleaners is recommended. High alkaline cleaning agents (>11.0pH) should not be used.

Note: alkaline cleaning agents may require neutralization following cleaning. Refer to the manufacturer's instructions to determine if neutralization is required, and follow the manufacturer's instructions for neutralization. It is critical alkaline cleaning agents be properly neutralized and rinsed from instruments.

Accessories: use of neutral detergents, enzymatic cleaners (7.0–8.0pH) soft bristled brushes, and soft pipe cleaners is recommended. Do not use alkaline cleaning agents as they cause significant damage and render accessories unusable.

The following solutions or tools should not be used on instruments or accessories:

- Saline solution.
- Solutions containing chlorine (e.g. bleach) or aldehydes (e.g. glutaraldehyde).
- Formalin, mercury, chlorides, bromides, iodides, or Ringer's solution.
- Metal brushes or scouring pads.

Cleaning and Rinsing Water

If available, softened tap water should be used. Critical water should be used for the final rinse step to prevent mineral deposits on surfaces (de-ionized water is recommended). One or more of the following processes may be used to treat water defined as critical water: ultra-filter (UF), reverse osmosis (RO), de-ionized (DI), or equivalent. Use of hard water should be avoided.

Use of Mineral Oil or Silicone-Based Lubricants

Mineral oil or silicone-based lubricants should not be used as they may not be removed by these cleaning instructions. These types of lubricants may coat microorganisms, prevent direct contact of steam with instrument surfaces, and hinder sterilization.

Processing Warnings/Cautions

- Do not allow soiled instruments or accessories to dry after use and prior to cleaning. Cleaning and subsequent sterilization may be hindered when blood or bloody solutions are allowed to dry on instruments or accessories.
- Medtronic reusable instruments are considered critical devices (i.e. critical contact) and must be thoroughly cleaned, then sterilized prior to initial use, or in adherence to these processing instructions before reuse.
- Soiled or used instruments should not be loaded into the accessories when cleaned in an automated washer-disinfector. Soiled instruments must be processed separately from accessories.
- Accessories can only be cleaned using enzymatic (neutral) cleaning agents.
- Universal precautions should be observed by health care facility personnel that work with contaminated or potentially contaminated instruments or accessories. Caution should be exercised when handling instruments with sharp points or cutting edges.
- Appropriate PPE should be worn when handling or working with contaminated or potentially contaminated materials, devices, and equipment. PPE includes gowns, masks, goggles, face shields, gloves, and shoe covers.
- Do not place heavy instruments on delicate instruments.
- Cool drafts from air ducts or other air currents should be avoided during the cooling phase to avoid post-sterilization moisture caused by rapid cooling.

Limitations on Processing

- Repeated processing has minimal effect on these instruments or accessories.
- End of life is determined by excessive wear and damage from normal use. Although treatment of instruments, accessories, materials used, and details of cleaning and sterilization have an important effect, there is no limit to the number of times instruments and accessories can be processed.
- See the EXAMINATION section of this document to determine if the device is at the end of its useful life.

INSTRUCTIONS - Instruments Only (Instructions for accessories are included after instructions for instruments)	
Bringing a new instrument into service	▪ Instruments are provided non-sterile and must be thoroughly cleaned, then sterilized before first use and any reuse. ▪ Remove packaging material prior to cleaning and sterilization.
Point of use	▪ Remove visible soil from instruments using non-shedding wipes. ▪ Instruments should be thoroughly cleaned within 30 minutes of use to minimize drying. ▪ If instruments cannot be processed immediately, keep instruments moist during transport. ▪ Place instruments in a tray of softened tap water, if available, or cover with damp towels.
Containment and transportation	▪ Immediately transport the tray containing covered instruments to a work area dedicated to further processing. ▪ Soiled instruments should be transported separately from non-contaminated instruments to avoid contamination of personnel and surroundings.
Preparation for cleaning	▪ Instruments that should be disassembled prior to cleaning and sterilization are provided with specific processing/disassembly/reassembly instructions. For additional details, refer to specific processing/disassembly/reassembly instructions. ▪ Remove occlusions prior to cleaning. Flush cannulated devices with tap or treated water to prevent drying of soil and/or debris on the inside. ▪ Instruments may be safely and effectively processed using manual or a combination manual and automated washer-disinfector (enzymatic or alkaline) cleaning procedures outlined in Tables 1 and 2.

Automated cleaning: enzymatic or alkaline cleaning agent	1. Rinse, flush, and scrub instruments under running (cold) tap water at a temperature of <43°C (110°F) for 30 seconds to 1 minute. Scrub instruments with appropriately-sized, soft bristle brushes to remove visible soil. Scrub inside any lumens or cavities and actuate instrument (if applicable). 2. Using tap water, prepare an enzymatic cleaning agent¹ (described after Table 1) or alkaline cleaning agent² (described after Table 2) according to the manufacturer's instructions, dilution recommendations, and temperatures in an appropriately-sized sonicator. 3. Place instruments in enzymatic cleaning agent, submerge completely, and ultrasonicate at 40–50kHz for a minimum of 15 minutes. Multiple ultrasonication cycles totaling a minimum of 15 minutes of ultrasonication time may be used. Ensure surfaces, lumens, and components of the device are in contact with cleaning agent. Note: ultrasonic cleaning is only effective if the surface to be cleaned is immersed in cleaning agent. Air pockets decrease efficacy of ultrasonic cleaning. Minimize air pockets or bubble formation by flushing lumens, cavities, crevices or springs with cleaning agent while the instrument is immersed in the ultrasonic tank. 4. Transfer instruments to the automated washer-disinfector and program the washer-disinfector with cycle parameters in Tables 1 (enzymatic) or 2 (alkaline). Ensure instruments that are cannulated and/or have blind holes are positioned in the washer-disinfector to allow for proper drainage, and pooling of wash water does not occur. Ensure cycle parameters in Tables 1 or 2 are properly programmed. 5. Upon unloading, visually inspect instruments, including lumens and cavities, to ensure visible soil is removed. If soil is still present, repeat the cleaning process or contact Medtronic immediately to arrange for disposal or replacement.
--	--

INSTRUCTIONS - Instruments Only
(Instructions for accessories are included after instructions for instruments)

Table 1: Enzymatic Cleaning Agent (7.0–8.9pH)

Cycle		Temperature set points	Minimum time (min:sec)	Concentration
Pre-wash		<43°C (110°F) Cold tap	2:00	N/A
Pulsed enzymatic wash ³	Wash	Hot tap	10:00	0.2–0.8mL/L
	Rinse	Hot tap	≥00:15	N/A
Detergent wash ⁴ (Neutral pH)		65°C (149°F)	5:00	0.2–0.8mL/L
Rinse		71°C (160°F)	15:00	N/A
Thermal (disinfection) Rinse (de-ionized)		≥90°C (194°F)	5:00	N/A
Dry		98.8°C (210°F)	≥30:00	N/A

Note: cycle validated using Prolystica™ Enzymatic Presoak and Cleaner¹ (1–4mL/L) for ultrasonication, Prolystica™ Ultra Concentrate Enzymatic Cleaner³, and Prolystica™ Ultra Concentrate Neutral Detergent⁴ for washer-disinfector use.

Table 2: Alkaline cleaning agent (8.8–11.8pH)

Cycle		Temperature set points	Minimum time (min:sec)	Concentration
Pre-wash		<43°C (110°F) Cold tap	2:00	N/A
Pulsed alkaline wash ²	Wash	Hot tap	10:00	2–6 mL/L
	Rinse	Hot tap	≥00:15	N/A
Alkaline detergent wash ²		60°C (140°F)	5:00	2–6 mL/L
Rinse		71°C (160°F)	15:00	N/A
Thermal (disinfection) Rinse (de-ionized)		≥90°C (194°F)	5:00	N/A
Dry		98.8°C (210°F)	≥30:00	N/A

Note: due to many variables involved with washer-disinfectors, health care facilities should properly install, calibrate, and verify the process (e.g. temperatures, times) used for their equipment. Washer-disinfector manufacturer recommendations should always be followed. When cleaning multiple devices in one cleaning cycle, ensure the manufacturer's maximum load is not exceeded. Cleaning validation was conducted using a full chamber load.

Disinfection: thermal	• Thermally decontaminate devices by a thermal final rinse within an automated washer-disinfector at ≥90°C (194°F) for 5 to 10 minutes. • Disinfectant solutions and chemicals are not required because the decontamination process includes cleaning followed by sterilization. • Thermal disinfection alone does not render instruments safe for patient use.
Drying: automated washer-disinfector	• Upon completion of washer-disinfector cycle, perform a visual inspection for dryness. • If moisture is observed on instruments after the washer-disinfector cycle, dry using clean, absorbent, non-shedding wipes and/or forced (medical grade) air, and/or dry heat oven at ≥118°C (245°F) for a minimum of 30 minutes. • For complex instruments (e.g. shaft-in sleeve, cannulation, and blind holes), a dry heat oven at ≥118°C (245°F) for a minimum of 30 minutes is required. Perform a visual inspection for dryness. • Instruments should be completely dry before packaged for sterilization.

Drying: manual	• Dry devices using clean, absorbent, non-shedding wipes, forced (medical grade) air, and/or a dry heat oven. Perform a visual inspection for dryness, and repeat if necessary. • Instruments should be completely dry before packaged for sterilization.
Inspection and reassembly	• If the instrument was disassembled prior to cleaning and sterilization, reassemble if applicable. • Carefully inspect the instrument for damage by performing a thorough examination as referred to the EXAMINATION section. • If the instrument is damaged, contact Medtronic immediately to arrange for disposal or replacement.
Packaging	• A single layer of Spunbond Meltblown Spunbond (SMS) polypropylene non-woven surgical sterilization wrap that is in conformance to ISO 11607-1 or equivalent must be used. • Individual instruments or dedicated sterilization trays must be double wrapped prior to sterilization. Wrap in accordance with local procedures per wrapping techniques listed in guidance document ANSI/AAMI ST79 or ISO/TS 16775. • Only instruments manufactured and/or distributed by Medtronic should be included. • Proceed to the Sterilization section.

INSTRUCTIONS - Accessories Only (Instructions for reusable instruments are provided in the previous table)

Bringing new accessories into service	• Accessories are provided non-sterile and must be thoroughly cleaned, then sterilized before first use and any reuse. • Remove packaging material prior to cleaning and sterilization.
Point of use	• Remove visible soil from accessories using non-shedding wipes. • Accessories should be thoroughly cleaned immediately after use to minimize drying.
Containment and transportation	• If accessories cannot be processed immediately, keep them moist during transport. • Use softened tap water on soiled accessories, if available, or cover with damp towels. • Immediately transport accessories to a work area dedicated to further processing.
Preparation for cleaning	• Accessories must be cleaned individually and not nested inside one another. Accessory lids must be removed and cleaned separately or cleaned in the open position if removal is not applicable. • Accessories may be safely and effectively processed using manual or combination manual and automated washer-disinfector (enzymatic only for accessories) cleaning procedures outlined in this document.
Automated cleaning: enzymatic	1. Rinse, flush, and scrub accessories under (cold) tap water <43°C (110°F) for a minimum of 30 seconds. 2. Transfer devices to the automated washer-disinfector and program the washer-disinfector with cycle parameters in Table 3. Ensure accessories are positioned in the washer-disinfector to allow proper drainage. Ensure inner portions of accessories are angled downward toward sprayer arms and cycle parameters in Table 3 are properly programmed. 3. Upon unloading, inspect devices, including any lumens and cavities/mated surfaces, to ensure visible soil is removed. If soil is still present, repeat the cleaning process or contact Medtronic to arrange for disposal or replacement.

Table 3: Enzymatic (neutral) cleaner (7.0-8.0pH)

Cycle		Temperature set points	Minimum time (min:sec)	Concentration
Pre-wash		<43°C (110°F) Cold tap	2:00	N/A
Pulsed enzymatic wash ⁵	Wash	60°C (140°F) Hot tap	4:00	0.2-0.8mL/L
	Rinse	60°C (140°F) Hot tap	≥00:15	N/A
Detergent wash ⁶ (Neutral pH)		60°C (140°F)	5:00	0.2-0.8mL/L
Rinse		71°C (160°F)	5:00	N/A
Thermal (disinfection) Rinse (de-ionized)		90°C (194°F)	5:00	N/A
Dry		98.8°C (210°F)	≥15:00	N/A

Note: cycle validated using Prolystica™ Ultra Concentrate Enzymatic Cleaner⁵ and Prolystica™ Ultra Concentrate Neutral Detergent⁶ for washer-disinfector use.

Note: due to many variables involved with washer-disinfectors, health care facilities should properly install, calibrate, and verify the process (e.g. temperatures, times) used for their equipment. Washer-disinfector manufacturer recommendations should always be followed. When cleaning multiple devices in one cleaning cycle, ensure the manufacturer's maximum load is not exceeded. Cleaning validation was conducted using a full chamber load.

Manual cleaning: enzymatic	1. Manually clean accessories when an automated washer-disinfector is not available. 2. Rinse, flush, and scrub accessories under cold tap water <43°C (110°F) for a minimum of 1 minute. 3. Scrub device with appropriately-sized, soft bristle brushes or soft pipe cleaners to remove visible soil. Scrub inside lumens or cavities and actuate devices, if applicable. 4. Scrub until visible soil is removed. 5. Using tap water, prepare an enzymatic cleaning agent⁷ according to manufacturer's instructions, dilution recommendations, and temperatures. Place devices in enzymatic cleaning agent, submerge completely, and soak for a minimum of 10 minutes. Ensure surfaces, lumens, and components are in contact with cleaning agent. 6. Remove devices from enzymatic cleaning agent and flush under running cold tap water. Lumens or cavities must be flushed in the water stream. Rinse with cold tap water for a minimum of 1 minute. 7. Using tap water, prepare a second enzymatic cleaning agent according to manufacturer's instructions, dilution recommendations, and temperatures in an appropriately-sized sonicator. 8. Place devices in enzymatic cleaning agent, submerge completely, and ultrasonicate at 40–50kHz for a minimum of 20 minutes. Multiple ultrasonication cycles totaling a minimum of 20 minutes of ultrasonication time may be used. Ensure surfaces, lumens, and components are in contact with cleaning agent. Note: ultrasonic cleaning is only effective if the surface to be cleaned is immersed in cleaning agent. Air pockets decrease efficacy of ultrasonic cleaning. Minimize air pockets or bubble formation by flushing lumens, cavities, crevices, or springs with cleaning agent while the device is immersed in the ultrasonic tank. 9. Remove devices from the sonicator and rinse using cold running tap water. Lumens or cavities must be flushed in the water stream. Rinse for a minimum of 10 minutes. Multiple cycles totaling a minimum of 10 minutes may be used. 10. Repeat rinsing as in step 9, this time with de-ionized water for an additional 10 minutes. Multiple cycles totaling a minimum of 10 minutes may be used. 11. Carefully inspect devices, including any lumens and cavities, to ensure visible soil is removed. If soil is still present, repeat the cleaning process or contact Medtronic immediately to arrange for disposal or replacement. Note: cycle validated using Prolystica™ Enzymatic Presoak and Cleaner⁷ (1–4mL/L).
Drying	Automated Washer-Disinfector ▪ Upon completion of washer-disinfector cycle, perform visual inspection for dryness. ▪ If moisture is observed on accessories, dry using clean, absorbent, non-shedding wipes and/or forced (medical grade) air until moisture is no longer present. Manual Dry devices using clean, absorbent, non-shedding wipes and/or forced (medical grade) air. Perform a visual inspection for dryness, and repeat if necessary.
Inspection and re-loading	▪ Inspect accessories for damage by performing a thorough examination as referred to the EXAMINATION section. ▪ If accessories or modules are damaged, contact Medtronic immediately to arrange for disposal or replacement. ▪ Load accessories with their applicable devices. ▪ Proceed to the Sterilization section.

STERILIZATION INSTRUCTIONS for both Instruments and Accessories

Sterilization

1. Double wrap individual instruments and/or dedicated accessories containing instruments with non-woven/moisture resistant surgical sterilization wrap.
2. Inspect packaging to ensure no rips, punctures, or seal failures are present in or on packaging prior to loading into the sterilizer.
3. Load instruments and accessories into the sterilizer by following sterilizer manufacturer's recommended loading procedures and load configurations.
4. Follow the sterilizer manufacturer's recommended procedures to program the sterilizer with any one of the sets of sterilization cycle parameters in Table 4.

Table 4: Sterilization cycle parameters

Cycle	Temperature set point	Minimum exposure time	Minimum dry time ⁸
Dynamic-air-removal	134°C (273°F)	3 Minutes	30 Minutes
Dynamic-air-removal	134°C (273°F)	18 Minutes	30 Minutes

⁸The minimum dry times were validated using sterilizers having vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer's recommendations.

Chamber size and chamber load differences may exist between industrial and health care facility sterilizer models. Sterilization parameters in Table 4 can be achieved in both health care facility and larger, industrial sterilizer models. Because of many variables involved in sterilization, health care facilities should calibrate and verify the sterilization process (e.g. temperatures, times) used for their equipment.

Note: steam for sterilization should be generated from water treated to remove total dissolved solids and non-condensable gases, filtered to remove contaminants and water droplets, and supplied via piping without dead legs or other stagnant zones where contamination may collect.

Caution: ethylene oxide (EO), gas plasma, gamma irradiation, chemical vapor, or dry heat sterilization methods are not recommended for sterilization of these devices. Steam/moist heat is the recommended method of sterilization.

Medtronic does not recommend or support Immediate-Use Steam Sterilization of instruments and accessories.

Additional information	Equipment for cleaning, disinfection, sterilization and drying must comply with ISO 17664, ISO 15883, and ISO 17665. Packaging materials must comply with ISO 11607. Instruments and accessories should be thoroughly cleaned, inspected, and sterilized before being returned to Medtronic. Any medical device that has been in contact with a patient known to have or suspected to have a prion or prion-related disease, such as Creutzfeldt-Jakob disease, should not be returned to Medtronic and should be quarantined and processed per guidelines for processing devices contaminated with high risk tissue at the health care facility. Contact Medtronic for information regarding proper disposal/destruction.
------------------------	---

Cleaning and sterilization information provided is in accordance with AAMI TIR12, AAMI TIR30, EN 285, ISO 17664, ISO 17665-1, ISO 15883-2, and ISO 15883-5. Processing instructions provided in this document were validated as being capable of preparing instruments for reuse. Sterilization validation demonstrated a sterility assurance level of 10^{-6} . Sterilization parameters were performed with accessories. Sterilization validation was conducted using an empty chamber load. A double wrapped, "worst case" perforated tray configuration was used in validation activities.

It remains the responsibility of the processor to ensure processing is performed using validated equipment and trained personnel in the processing facility to achieve desired results. This normally requires validation and routine monitoring of the process. Deviation by the processor from instructions provided should be properly evaluated for effectiveness and potential adverse consequences. Some devices have specific disassembly/assembly instructions. In this case, refer to instructions provided with the instrument for additional instructions. Users must establish appropriate processing protocols for reusable medical devices using recommendations of the device manufacturer and cleaning agent manufacturer.

STORAGE AND SHELF LIFE

For non-sterile instruments and accessories, store in dry, clean conditions at ambient room temperature. Non-sterile instruments and accessories are reusable and do not have a specified shelf life. End of life is determined by excessive wear and damage from normal use. Refer to the Limitations on Processing section for more information.

For sterilized instruments and accessories, ensure instruments and accessories are dry before storing. Store in clean, dry conditions at ambient room temperature. Shelf life of sterilized instruments and accessories is dependent upon the legally marketed sterilization wrap. The health care facility should establish a shelf life for sterilized instrumentation based on the type of sterilization wrap used and recommendations of the sterilization wrap manufacturer.

PRODUCT COMPLAINTS

To report product problems, contact Medtronic.

Patients in the European Union experiencing a serious incident in relation to the device should contact Medtronic and the competent authority of the Member State in which they are established.

MRI INFORMATION

Medtronic instruments and accessories are not intended to be used in the magnetic resonance (MR) environment. As such, Medtronic instruments and accessories were not evaluated for safety and compatibility in the MR environment. Therefore, safety of Medtronic instruments and accessories in the MR environment is unknown.

FURTHER INFORMATION

Recommended directions for use (surgical operative techniques) are available at no charge upon request. If further information is required, contact Medtronic.

©2020 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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AUSTRALIAN SPONSOR:
Medtronic Australasia Pty Ltd
5 Alma Road
Macquarie Park, NSW 2113
Australia

Medtronic Sofamor Danek USA, Inc.
1800 Pyramid Place
Memphis, TN 38132
Telephone: 800 933 2635 (USA)
901 396 3133 (Outside USA)
Fax: 901 396 0356

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel: + 31 45 566 80 00

EXPLANATION OF SYMBOLS

Authorized representative in the European Community
Упълномощен представител в Европейската общност
Autorizovaný zástupce pro Evropské společenství
Autoriseret representant i Det Europæiske Fællesskab
Autorisierter Repräsentant in der Europäischen Gemeinschaft
Εξουσιοδοτημένος αντιπρόσωπος στην Ευρωπαϊκή Κοινότητα
Representante autorizado en la Comunidad Europea
Volikud esindeja Euroopa Ühendusea Valtuudetü edustaja Euroopan yhteisön alueella
Représentant autorisé dans la Communauté européenne
Ovlašteni predstavnik u Evropskoj zajednici
Hivatalos képviselő az Európai Közösségekben
Rappresentante autorizzato nella Comunità europea

The device complies with European Directive MDD 93/42/EEC
Устройство соответствует на Европейска Директива 93/42/ЕО относно медицинските изделия
Zařízení splňuje požadavky evropské směrnice MDD 93/42/EHS
Enheden overholder Det Europæiske Direktiv MDD 93/42/EF
Die Vorrichtung erfüllt die Anforderungen der europäischen Richtlinie 93/42/EWG über Medizinprodukte
Το προϊόν προϋπὸν συμμορφώνεται με την Ευρωπαϊκή Οδηγία περί ιατροτεχνολογικών προϊόντων 93/42/ΕΟΚ
El dispositivo cumple con la Directiva europea 93/42/CEE relativa a los productos sanitarios
Seade vastab Euroopa meditsiiniseadmeid käsitlevale direktiivile 93/42/EMO
Tämä laite noudattaa eurooppalaista lääketieteellisedirektiivin 93/42/ETY
Le dispositif est conforme à la directive européenne MDD 93/42/CEE
Proizvod je usklađen s evropskom Direktivom MDD 93/42/EEZ
Az eszköz megfelel az orvostechnikai eszközökről szóló 93/42/EGK európai irányelvnek
Questo dispositivo è conforme alla Direttiva europea MDD 93/42/CEE

Igalizotais atstovas Europos Bendrijoje
Pārtvarotais pārstāvis Eiropas Kopienā
Geautoriseerd vertegenwoordiger in de Europese Gemeenschap
Autoriseret representant i Det europeiske fællesskab
Autoryzowany przedstawiciel we Wspólnocie Europejskiej
Representante autorizado na Comunidade Europeia
Reprezentantă autorizată în Comunitatea Europeană
Уполномоченный представитель в Европейском сообществе
Autorizovaný zástupce pro Evropské společenstvo
Pooblaščený predstavnik v Evropski skupnosti
Ovlašteni predstavnik u Evropskoj zajednici
Auktoriserad representant inom EG
Avrupa Topluluğundaki yetkili temsilci

Šis prietais atitinka Europos Direktyvą 93/42/EEB dėl medicinos prietaisų
Ierice atbilst Eiropas Padomes Medicīnas ierīču direktīvas 93/42/EEK prasībām
Dit product voldoet aan de Europese Richtlijn MDD 93/42/EEG
Enheten er i samsvar med EU-direktiv MDD 93/42/EF
To urządzenie jest zgodne z dyrektywą europejską MDD 93/42/EWG
O dispositivo está em conformidade com a Diretiva Europeia MDD 93/42/CEE relativa a dispositivos médicos.
Dispozitivul este în conformitate cu Directiva europeană 93/42/CEE privind dispozitivele medicale
Устройство соответствует Европейской директиве MDD 93/42/EEC
Pomůcka splňuje ustanovenia směrnice EÚ o zdravotnických pomůckách 93/42/EHS
Připomoček je v skladu z evropsko Direktivo (93/42/EGS) o medicinskih pripomočkih
Ovo sredstvo je usklađeno sa Evropskom direktivom za medicinska sredstva 93/42/EEZ
Detta instrument uppfyller Europarådets direktiv MDD 93/42/EEG
Bu cihaz MDD 93/42/EEC sayılı Avrupa Direktifi ile uyumludur

The device complies with European Directive MDD 93/42/EEC
Устройството съответства на Европейска Директива 93/42/EО относно медицинските изделия
Zařízení splňuje požadavky evropské směrnice MDD 93/42/EHS
Enheden overholder Det Europæiske Direktiv MDD 93/42/EF
Die Vorrichtung erfüllt die Anforderungen der europäischen Richtlinie 93/42/EWG über Medizinprodukte
Το παρόν προϊόν συμμορφώνεται με την Ευρωπαϊκή Οδηγία περί ιατροτεχνολογικών προϊόντων 93/42/ΕΟΚ
El dispositivo cumple con la Directiva europea 93/42/CEE relativa a los productos sanitarios
Seade vastab Euroopa meditsiiniseadmeid käsitlevale direktiivile 93/42/EMO
Täna laite nõustab euroopaalaista laakintalaitedirektiivis 93/42/ETY
Le dispositif est conforme à la directive européenne MDD 93/42/CEE
Proizvod je usklađen s europskom Direktivom MDD 93/42/EEZ
Az eszköz megfelel az orvostechnikai eszközökről szóló 93/42/EGK európai irányelvnek
Questo dispositivo è conforme alla Direttiva europea MDD 93/42/CEE

Šis prietaisas atitinka Europos Direktyvą 93/42/EEB dėl medicinos prietaisų lerice atbilst Eiropas Padomes Medicinas ierību direktīvas 93/42/EEK prasībām
Dit product voldoet aan de Europese Richtlijn MDD 93/42/EEG
Enheden er i samsvær med EU-direktiv MDD 93/42/EF
To urządzenie jest zgodne z dyrektywą europejską MDD 93/42/EWG
O dispositivo está em conformidade com a Diretiva Europeia MDD 93/42/CEE relativa a dispositivos médicos.
Dispozitivul este în conformitate cu Directiva europeană 93/42/CEE privind dispozitivele medicale
Устройство соответствует Европейской директиве MDD 93/42/EEC
Pomůcka splňuje ustanovenia směrnice EÚ o zdravotnických pomůckách 93/42/EHS
Pripomoček je v skladu z evropsko Direktivo (93/42/EGS) o medicinskih pripomočkih
Ovo sredstvo je usklađeno sa Evropskom direktivom za medicinska sredstva 93/42/EEZ
Delta instrumenti uppfyller Europarådets direktiv MDD 93/42/EEG
Bu cihaz MDD 93/42/EEC sayılı Avrupa Direktifi ile uyumludur

METRx® MICROSCOPE

0381036 Rev. B

2016-10-26

Medtronic

* 0 3 8 1 0 3 6 - B *

ENGLISH

IMPORTANT INFORMATION ON THE METRx® MICROSCOPE SYSTEM

DESCRIPTION

The METRx® Microscope System is composed of a microscope, camera, an integrated video system, a video monitor and recorder, and various cannula, dilators, extension lenses, and associated instruments.

No warranties, express or implied are made. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

INDICATIONS

The METRx® Microscope is indicated for visualization of the surgical field in any area of the body cut open during a surgical procedure. When used in the cervical, thoracic, or lumbar spine either from an anterior or posterior direction, for example, the METRx® Microscope and accessories are intended to aid the surgeon's visualization of the surgical area and allow him/her to perform any type of surgical spinal procedure such as herniated disc repair, visualization of the circumferential decompression of the nerve roots, aiding in the search and removal of nucleus material, spinal fusion, or insertion of spinal implants. Other examples of generic surgical use of the METRx® Microscope would be for use in the knee, ankle, shoulder, hand, wrist, and temporomandibular joint (TMJ).

CONTRAINDICATIONS

The METRx® Microscope and accessories have no known contraindications intrinsic to the device. No part of the microscope itself or its accessories should ever be used in a cutting or tearing action, i.e., never use the Microscope as an instrument. The optical device or its extension lenses should not be used to provide access to the surgical field. The device should not be inserted into body cavities, hollow organs, or natural body openings.

There are no other known risks associated with the use of the device outside of the normal and expected risks of surgery. The microscope should not be used non-sterile or in the presence of an infectious disease process.

DIRECTIONS FOR USE

Specific instructions for use depend on patient considerations. Therefore, Medtronic cannot provide a surgical procedure applicable to all situations. Any available surgical procedure brochure or manual for the METRx® System Microscope device and accessories should be reviewed prior to use. The only critical directions for use are to insert the cannula or dilators and position them in the surgical wound prior to insertion of the microscope extension lenses and inner field visualization of the surgical site. Never use the microscope extension lenses to provide access or as a surgical instrument in the surgical field. Once visualization assistance is obtained, the surgeon can complete the planned surgical procedure.

POTENTIAL ADVERSE EFFECTS

Risks possibly associated with the use of the METRx® Microscope and accessories are similar to those associated with any surgery to the planned area of instrument use. The most frequently stated risks are bleeding, neurological damage, damage to the surrounding soft tissue, and infection. Each of these risks has also been used to describe the risks associated with conventional surgical intervention. Additional risks associated with the use of METRx® Microscope, other than those described for spinal surgery in general, may be instrument malfunction, such as bending, fragmentation, loosening, and/or breakage (whole or partial). Breakage of the tip in the patient may increase surgical time, since this instrument should not be implanted. Also, the surgery may not be effective. Similar risks are associated with the system use in other parts of the body.

Additional risks are attendant to surgery and the use of anesthesia, etc., and are not directly related to the use of the microscope and accessories. These include, but are not limited to, pneumonia, phlebitis, embolism, wound infection, and blood loss with or without anemia.

WARNINGS AND PRECAUTIONS

A successful result is not always achieved in every surgical case. This fact is especially true in orthopaedic or neuro-surgery cases where many extenuating circumstances may compromise the results.

In the event of technical complications, the surgical technique can be converted to an open procedure and the surgery completed.

Preoperative and operating procedures, including knowledge of surgical techniques are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection of the patient and the compliance of the patient will greatly affect the results.

In addition, the following should be considered:

1. This device is a delicate instrument. It should NOT be dropped, bent at a sharp angle, or exposed to any type of gamma radiation. Tip fracture fragmentation and optic damage may result if the microscope is not handled carefully.
2. Additional microscopes and accessories should be available at the time of surgery in case of possible contamination due to mishandling or removing the devices from the sterile field.
3. Components of the system should be thoroughly inspected during cleaning prior to surgery for possible damage.
4. Proper, secure component connections must be made to ensure proper functioning of the optical, irrigation, and aspiration aspects of the device.

Caution – high temperature.

Typically, a light source is used in conjunction with the METRx® System. To reduce the risk of patient injury or damage to materials that may come into contact with this illumination system, always ensure the heat shield is in place covering the metal cable connectors. It is recommended that this illumination system not be used with light sources rated over 300 Watts or with fiber optic light cables larger than 5mm. However, due to inherent differences in light and power output among light sources, we recommend the connector be monitored for elevated temperatures. Use caution when touching the metal connectors when disconnecting the light cable after use.

CLEANING PROCEDURE

Exterior cleaning of the microscope and camera head is essential prior to any sterilization procedure. The camera head should be removed from the microscope prior to cleaning. To remove blood, organic matter, and irrigation solutions, all surfaces of the microscope and camera head should be cleaned with a mild detergent and water. Rinse with distilled, demineralized, or pyrogen-free water. Dry thoroughly. Without the removal of all contaminants from the surface, the sterilization medium will not contact the surfaces.

Warning: do not use ultrasonic cleaner or abrasives during the cleaning process.

STERILIZATION:

The microscope, light cable, non-sterile instruments, and instruments which are re-usable are recommended to be steam sterilized by the hospital using one of the following methods:

Some cameras (9580500, 9580501) are autoclavable following the same recommendations. If the products described in this document are sterilized by the hospital in a tray or case, they must be sterilized in a tray or case provided by Medtronic.

Some accessories and instruments are supplied sterile and non-reusable. Sterile product will be clearly labeled as such on the package label. The sterility of the product supplied sterile can only be assured if the packaging is intact.

Note: the following note applies to the process parameter identified with the * below. For use of this product and instruments outside the US, some non-US health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical instruments that could come into contact with the central nervous system.

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME
Steam	Pre-Vacuum	270°F (132°C)	4 Minutes
Steam	Gravity	250°F (121°C)	60 Minutes
Steam*	Pre-Vacuum*	273°F (134°C)*	20 Minutes
Steam*	Gravity*	273°F (134°C)*	20 Minutes*

The microscope, camera, and light cable can also be EtO sterilized. The following EtO cycle is recommended:

Preconditioning parameters:

Temperature:	55 ± 2° C
Relative Humidity:	≥ 35%
Vacuum:	21 ± 1 in Hg (508-550mm Hg)
Preconditioning time:	1 Hour

Sterilization parameters:

Ethylene Oxide Carrier:	Oxyfume 2002
Temperature:	55 ± 2° C
Relative Humidity:	≥ 35%
Pressure:	19 ± 1 PSIG (2.25-2.39 bars)
Ethylene Oxide Concentration:	736 mg/L
Gas Exposure Time (Full Cycle):	4 Hours
Aeration:	11 Hours at 54°C minimum

The microscope and camera should be thoroughly cleaned prior to sterilization.

The integrated video system, video monitor, and recorder are also reusable and supplied non-sterile. These components should not be placed in the surgical field.

Caution: microscopes and camera cannot be steam autoclaved unless "autoclavable" is engraved on the microscope or camera body. This method would otherwise permanently damage the optical components.

Caution: do not immerse or rinse instruments in cold water or any other liquid to accelerate cooling.

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MRI INFORMATION

Medtronic instruments are not intended to be used in the magnetic resonance (MR) environment. As such, Medtronic instruments have not been evaluated for safety and compatibility in the MR environment. Therefore, the safety of the Medtronic instruments in the MR environment is unknown.

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© 2016 Medtronic Sofamor Danek, Inc. All rights reserved.

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* 0 3 8 1 2 4 2 - C *

ENGLISH

IMPORTANT INFORMATION ON THE MAST QUADRANT® RETRACTOR SYSTEM

DESCRIPTION

The MAST QUADRANT® Retractor System is a tubular-based retraction system, designed to provide surgeons with the freedom to retract tissue through any combination of distracting or articulating the blades. The MAST QUADRANT® Retractor System includes instruments used to access the spine by dilating the overlying tissues, as well serving as a retracting device to maintain the access. The system may be used in conjunction with microscopes, light sources, cameras, or other visualization aids. No warranties, express or implied, are made. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

INDICATIONS

The MAST QUADRANT® Retractor System is intended to provide surgeons with instruments such as dilators, retractors, light sources and pedicle access needles used to perform a variety of spinal fixation procedures utilizing a minimally invasive approach.

The MAST QUADRANT® Retractor System is indicated for visualization of the surgical field in any area of the body cut open during a surgical procedure. When used in the cervical, thoracic, or lumbar spine either from an anterior or posterior direction, for example, the MAST QUADRANT® Retractors and accessories are intended to aid the surgeon's visualization of the surgical area and allow him/her to perform any type of surgical spinal procedure such as herniated disc repair, visualization of the circumferential decompression of the nerve roots, aiding in the search and removal of nucleus material, spinal fusion, or insertion of spinal implants.

CONTRAINDICATIONS

The MAST QUADRANT® Retractor System has no known contraindications intrinsic to the device. No part of the device should ever be used in a cutting or tearing action. The device should not be inserted into body cavities, hollow organs, or natural body openings. There are no other known risks associated with the use of the device outside of the normal and expected risks of surgery.

DIRECTIONS FOR USE

Specific instructions for use depend on patient considerations. Therefore, Medtronic cannot provide a surgical procedure applicable to all situations. Any available surgical procedure brochure or manual for the MAST QUADRANT® Retractor System and for all instruments should be reviewed prior to use. The only critical directions for use are to insert the cannula or dilators and position them in the surgical wound prior to insertion of the device. Once visualization assistance is obtained, the surgeon can complete the planned surgical procedure.

POTENTIAL ADVERSE EFFECTS

Risks possibly associated with the use of the MAST QUADRANT® Retractor System are similar to those associated with any surgery to the planned area of Retractor use. The most frequently stated risks are bleeding, neurological damage, damage to the surrounding soft tissue, and infection. Each of these risks has also been used to describe the risks associated with conventional surgical intervention. Additional risks associated with the use of the MAST QUADRANT® Retractors, other than those described for spinal surgery in general, may be instrument malfunction, such as bending, fragmentation, loosening, and/or breakage (whole or partial). Also, the surgery may not be effective. Similar risks are associated with the system use in other parts of the body.

Additional risks are attendant to surgery and the use of anesthesia, etc., and are not directly related to the use of the instruments. These include, but are not limited to, pneumonia, phlebitis, embolism, wound infection, and blood loss with or without anemia.

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1. This device is a delicate instrument. It should not be dropped, bent at a sharp angle, or exposed to any type of gamma radiation.
2. Additional accessories should be available at the time of surgery in case of possible contamination due to mishandling or removing the devices from the sterile field.
3. Components of the system should be thoroughly inspected during cleaning prior to surgery for possible damage.
4. Proper, secure component connections must be made to ensure proper functioning of the optical, irrigation, and aspiration aspects of the device.

Caution: high temperature.

Typically, a light source is used in conjunction with the MAST QUADRANT® Retractors.

To reduce the risk of patient injury or damage to materials that may come into contact with this illumination system, always ensure the heat shield is in place covering the metal cable connectors. It is recommended that this illumination system not be used with light sources rated over 300 Watts or with fiber optic light cables larger than 5mm. However, due to inherent differences in light and power output among light sources, we recommend the connector be monitored for elevated temperatures. Use caution when touching the metal connectors when disconnecting the light cable after use.

For US audiences only

Caution: federal law (USA) restricts these devices to sale by or on the order of a physician.

CLEANING PROCEDURE

Exterior cleaning is essential prior to any sterilization procedure. Dry thoroughly. Without the removal of all contaminants from the surface, the sterilization medium will not contact the surfaces.

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STERILIZATION

The non-sterile instruments and instruments which are re-usable are recommended to be steam sterilized by the hospital using one of the following methods:

If the products described in this document are sterilized by the hospital in a tray or case, it must be sterilized in a tray or case provided by Medtronic.

Some accessories and Retractors are supplied sterile and non-reusable. Sterile product will be clearly labeled as such on the package label. The sterility of the product supplied sterile can only be assured if the packaging is intact.

Note: the following note applies to the process parameter identified with the * below. For use of this product and Retractors outside the US, some non-US health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical Retractors that could come into contact with the central nervous system.

Method:	Steam	Steam	Steam
Cycle:	Pre-Vacuum	Gravity	Gravity
Temperature	270°F (132°C)	250°F (121°C)	273°F (134°C)*
Exposure Time:	4 minutes	60 minutes	20 minutes*

The MAST QUADRANT® Retractors should be thoroughly cleaned prior to sterilization.

Caution: do not immerse or rinse light sources used with these instruments in cold water or any other liquid to accelerate cooling.

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INSTRUCTIONS FOR USE - REUSABLE SURGICAL INSTRUMENTS

EN

IMPORTANT INFORMATION FOR REUSABLE SURGICAL INSTRUMENTS FROM GAUTHIER BIOMEDICAL

SCOPE

This instruction for use will apply to reusable manual surgical instruments that are torque handles, ratcheting handles and ratcheting torque handles.

INTENDED USE

This instrument is intended for use in surgical procedures on bone and tissue by a physician trained in its intended use. Some instruments may be used in conjunction with other instruments.

This instrument may be connected to shaft type instruments such as awls, taps, screw drivers and similar instruments used to prepare for installation, driving or adjusting fasteners or constructs into bone with the appropriate instrument connector.

Device is sold non-sterile. Sterilize before use per manufacturer's reprocessing instructions below. This instrument may be used in orthopedic surgical procedures which are also known as bone surgical procedures in Japan.

Attachment Mechanism

Actuate attachment mechanism (or release collar if present) at front of handle, then insert shaft type instrument into connector while rotating slightly until it clicks into place. Check that shaft type instrument assembles readily with the mating connector on the torque handle.

Pull on shaft type instrument to ensure that it is securely connected.

Push on release collar to release shaft type instrument from the torque handle.

Mating Geometry

Instruments intended to mate with other instruments or the implant during use should be verified prior to use. Any mating geometry present on the instrument should be verified to fit with the appropriate mating parts before use.

Torque Mechanism

The torque handle is intended to limit torque transmission to a preset torque value or to indicate torque.

Torque handles may be bi-directional or uni-directional depending on the mechanism. Verify direction of torque mechanism before use.

Verify torque value before use. Ensure that a backup torque handle is available should the main torque handle be out of calibration.

Calibration: Device must be sent to manufacturer for evaluation or calibration every 6 months, unless otherwise specified on the device labeling. Device must be sent to manufacturer for evaluation or calibration when torque handle is determined to be out of calibration.

Ratcheting Mechanism

For Ratcheting Handle: handles with a ratcheting mechanism may have a DIRECTION CONTROL CAP with which to select the ratcheting direction or locked position. It is important for the user to familiarize themselves with the device functions before use in the surgery.

If a direction control cap is present, while holding the handle, rotate the direction control cap through the various detent positions. There could be two or three detents.

Each detent position will be one of the following:

- (1) Drive clockwise and ratchet counterclockwise
- (2) Locked in both directions or
- (3) Drive counter clockwise and ratchet clockwise

If no direction control cap is present, the instrument will drive in one direction and ratchet in the other direction.

It is important to verify the ratcheting handle drive direction before use.

DO NOT IMPLANT THE INSTRUMENTS.

Any attempt to change or modify this instrument will void the warranty.

WARNINGS

If the instrument slips, is dropped, mishandled or used improperly the patient and health care personnel may be injured by sharp edges.

Failure to follow proper maintenance, handling, cleaning, sterilization and storage procedures could make the instrument unfit for use and even be dangerous to the patient or health care personnel.

Impacts applied to the instrument not intended for impact or excessively high impact forces could cause the instrument to become damaged.

For Torque Handle with An Impact Cap: If an Impact cap is present on the handle, the impact cap will be engraved with the words "IMPACT CAP". If the words "IMPACT CAP" are not engraved on a cap then it is not an impact cap. Handles often have caps that are not impact caps.

If an impact cap is present on the instrument then the instrument may be impacted with light to medium impacts and only on the impact cap. Impact may only be applied with a surgical mallet

or similar instrument.

A medium impact is defined as an impact applied with a mallet not exceeding 1.65lb (748 gms) from a distance of no more than (14) inches (356 mms) with a medium speed. Medium speed is defined as approx. ½ - 1 sec from start to end of impact stroke.

A light impact is defined as an impact applied with a mallet not exceeding 1.65lb (748 gms) from a distance of no more than 4 inches (102 mms) with a medium speed. Medium speed is defined as approx. ½ - 1 sec from start to end of impact stroke.

Sharp edges could cut or injure health care personnel and patient. Chips, fragments or components could dislodge or break off the instrument and become implanted in the patient. The instrument is not intended to be implanted.

Rod bending or rod cutting instruments have the specific risks of producing broken bits of material. These instruments could apply very high forces and sometimes bits of material can be produced when rods or other materials are bent or cut. This is why it is hazardous to use them in situ. Any material left in the body of the patient could cause an allergic reaction or an infection. The physician and rest of healthcare personnel should wear PPE (personal protection equipment) like mask, eye protection (goggles), gloves.

The physician should take extreme care when operating near vital organs, nerves or blood vessels as slippage, excessive force or a combination of the two could cause harm.

PRECAUTIONS

Take care while applying high forces as this could dislodge or damage (implants/constructs) (devices, particularly hooks).

Do not expose the instrument to temperatures exceeding 135°C as this could damage the instrument.

Take care to make sure that instruments are well maintained and remain fully functional.

Handle instruments with care. Store the instrument in a clean, dry space that avoids direct sunlight and extreme temperatures.

Do not redesign or modify the instruments in any way as this could damage the instrument, make it dangerous to operate and it may not function as intended.

Do not use the instrument for any use other than that which it is intended for.

Periodically review instrument function and make use of repair and replacement services.

Do not use damaged instruments or instruments that are not fully functional. Additional back-up instruments should be available.

The physician must be knowledgeable of pre- and post-operative procedures in order to properly use the instrument.

Proper patient selection and their adherence to recommended procedures will influence outcomes.

Care should be taken while operating on pediatric patients since they are more vulnerable to stresses from the manipulation of bone and tissue.

If the instrument has a measuring function ensure that all necessary markings are clearly visible.

POSSIBLE ADVERSE EFFECTS

Damage to nerves, soft tissue, organs or joints could result in paralysis and pain.

The health care personnel or patient could become infected from improperly cleaned or sterilized instruments. The patient could also become infected due to surgical complications. The patient could experience pain, discomfort or abnormal sensations when the instrument is present.

Nerves are fragile and could be damaged during surgery.

A tear in the membrane that surrounds the brain and spinal cord known as Dura can leak as an unintentional consequence of surgery.

Blood vessels, nerves and organs could be struck and damaged by an accidentally dropped or misplaced instrument.

Injury may be caused by sudden unintentional release of clamped or spring loaded instrument components.

Skin or gloves of health care personnel could be cut by sharp edges.

Bone fracture especially in the case of a weak or deformed spine could occur.

Injury to patient or healthcare personnel may result from the instrument components breaking apart, becoming disassembled.

The health care personnel, especially the physicians' training and expertise in the surgical procedure, determine the usage of the instrument. The patient outcomes are not always positive and some surgeries are not successful. A positive patient outcome (successful surgery) is dependent on a variety of patient factors.

Severe inflammatory reactions from endotoxin contamination of surgical instruments may result.

Used ultrasonic cleaner liquid may be contaminated by bacteria or the cleaning fluid could result in endotoxin contamination.

Alkaline cleaners are corrosive.

Alcohols are not recommended for sterilizing medical materials because they do not have the ability to effectively kill spores and they cannot penetrate protein-rich materials.

Alcohols tend to swell and harden rubber and some plastic tubings after prolonged and repeated use. Alcohols also tend to bleach rubber and plastic tiles.

Hinged instruments must be fully opened when loaded in an automated washer to allow adequate contact with the detergent solution.

Do not stack instruments in the automated washer as this will reduce surface area exposed to contact with the detergent solutions, flushing and rinsing liquids.

Instruments should be disassembled when possible to improve cleaning.

PHYSICIAN NOTE: Since the physician is the person between the company and the patient, the important medical information in this document should be communicated to the patient.

CAUTION: FEDERAL (U.S.) LAW RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN ONLY.

This instrument should only be used by physicians that know and understand its intended use, the mating instruments and the surgical procedures that this instrument is to be used for.

PACKAGING

Check to make sure the packaging is not open nor damaged when received. Check product to make sure there is no damage and product is in working order prior to use.

Damaged packages should be returned to local distributor.

Remove instrument from packaging then clean and sterilize instrument prior to use in surgery. Instruments must be cleaned and sterilized prior to returning product to distributor. Health care personnel or the responsible person should include a note to notify the recipient that the instrument has been cleaned and sterilized and therefore is not a biohazard.

Instrument cases must be wrapped with a sterilization wrap to create a sterile barrier and maintain sterility of instruments within the case.

MAINTENANCE

Carefully examine instruments visually and functionally to ensure they are fully functional and safe before using them in surgery.

Lubricate instruments in order to maintain the action of moving parts and to protect the instrument from staining and rusting during reprocessing and storage, follow instructions in the Lubrication section of these instructions.

Store the instrument in a clean, dry space that avoids direct sunlight and extreme temperatures.

EXAMINATION

The user must examine instruments both visually and functionally before using them in surgery. Examine the instrument to make sure that the following features are visually and functionally acceptable and that none of the components are missing:

working surfaces

pivots

racks

spring when accessible or viewable

torsional components

cleanliness of holes or lumens

presence of cracks

bent components or entire instrument

deformed

distorted

Do not use instruments that appear too worn to use, damaged, missing components or otherwise non-functional.

Visual Inspection

Check to ensure:

All markings are clearly visible.

Instrument does not have any cracks.

There are no stains, discoloration or rust.

Check to ensure that handle to shaft connection is secure.

There are no cuts, nicks, tears, gouges or other damage to the silicone.

Insulation is not damaged.

The (working end) functional features are not damaged. Sharp edges should be sharp.

Threads appear to be free of damage. Check fit with mating part.

Ensure no parts are missing, loose or damaged. These parts include but are not limited to springs, screws and pins.

No damage that would inhibit or affect mating function. Check fit with mating part.

Check to ensure lumens are intact and free of obstructions.

Functional Inspection

Check to ensure: (where present)

Moving parts move freely (smoothly without being sloppy).

Springs function as intended to return the lever, handle or component it is connected to back to its original position.

Retention tabs are intact and fully functional.

Instrument mates and functions with compatible mating parts.

Ball detents are intact and function to hold parts together.

Sharp edges appear intact and sharp.

Tips close together properly.

Visible ratchet teeth are not damaged or missing. Verify ratcheting mechanism is fully functional by ratcheting the instrument at least one full rotation. When a direction control cap is present, verify driving and ratcheting directions in all positions. Mating component may be required to verify ratcheting and driving functions.

Driver tips are not excessively worn. Check mating with appropriate part.

Do not bend or cut rods in situ or in vivo as fragments could break off into the patient or rod bender could slip or accidentally spring open and strike implants causing harm to patient.

Bending, unbending or re-bending could induce fatigue and possible failure in the rod implant.

Hand-bendable rod templates can be used to determine ideal rod implant shape and positioning (placement) prior to bending and placing the rod implant.

There is no exact limit to the life of the instrument. The instrument may be used until it fails visual or functional inspection.

REPROCESSING – GENERAL INSTRUCTIONS (CONSIDERATIONS)

Cleaning is defined as the removal of organic material or soil, microorganisms and lubricants.

Proper cleaning will accomplish the following:

Reduce to a minimum or eliminate the transfer of soil from one patient to another.

Helps ensure that the following steps of reprocessing are successful. Reduces the level of soil enough to help make sterilization successful and therefore helps the entire reprocessing process be successful.

Sterilization should promptly follow the cleaning steps in order to accomplish successful reprocessing.

NOTE: Instruments requiring disassembly should be taken apart prior to cleaning.

Blood borne Pathogens

Universal precautions for handling this device after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030, Occupational Exposure to Blood borne Pathogens.

Instruments used on patients with suspected or confirmed infectious disease especially Creutzfeldt-Jakob Disease (CJD) or Transmissible Spongiform Encephalopathy (TSE) infected tissue must be thoroughly sterilized according to local laws and procedures.

These instructions do not include recommendations for the inactivation of TSE agents also known as prions. The user may refer to the World Health Organization (WHO) Infection Control Guidelines for TSE for instructions on the inactivation of prions.

Automated Cleaning

For thorough cleaning it is recommended that automated cleaning be preceded by manual cleaning steps as automated cleaning alone may not be enough.

DISINFECTION

Thermal disinfection may be performed as an additional step to process devices, but devices must be manually cleaned and terminally sterilized prior to their use in surgery.

Manual Cleaning

It is recommended that enzymatic cleaners be used to remove organic soil. Disinfection solutions may not have the cleaning properties of an enzymatic cleaner.

These cleaning instructions were validated using an enzymatic cleaner.

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The following cleaning agents, solutions, or tools should NOT be used:

- Saline solution
- Alkaline cleaning agents
- Solutions containing chlorine (e.g., bleach) or aldehydes (e.g., glutaraldehyde)
- Formalin, mercury, chlorides, bromides, iodides, or ringers solution
- Metal brushes or scouring pads

The following cleaning agents and tools are recommended for use:

- Soft bristled brushes and pipe cleaners of appropriate sizes
- Neutral pH enzymatic cleaners
- Cloths and towels used for cleaning must be clean and lint free.
- Use only low foaming cleaning agents so that instruments are visible in the cleaning solution.

Discard worn brushes as these will not clean effectively.

Use appropriate size brush for lumen being cleaned because if it is too small it will not clean thoroughly and if it is too big it will not fit into the lumen.

Brushes, towels and other cleaning tools should themselves be cleaned and disinfected or sterilized so that they do not add bioburden to the instruments .

If using highly alkaline cleaners the instrument must be properly neutralized and then thoroughly flushed and rinsed with de-ionized water.

Use of alkaline chemicals could cause rusting and pitting.

Highly aggressive agents such as Sodium hydroxide (NaOH) or Sodium hypochlorite (NaOCl) commonly known as bleach may discolor or corrode instruments.

Use of cleaners with rust inhibitors could help prevent or reduce the formation of rust. Follow manufacturer's directions for temperature, dilution and application.

Sort instruments by similar metal for subsequent processing so that electrolytic (galvanic corrosion) due to contact between dissimilar metals will not occur.

Plated instruments should not be cleaned in an ultrasonic cleaner since the ultrasonic vibration and the presence of other sharp instruments may crack or rupture the plating.

Rinse cleaning agents thoroughly off the instruments to prevent accumulation of residue.

Note: features that are difficult to clean include but are not limited to blind holes, lumens, rough surfaces, teeth, springs, ball detents, and mated surfaces with narrow gaps.

Cleaning and Rinsing Water

Use softened tap water when available. The use of de-ionized water for the final rinse step will prevent mineral deposits or water stains being left behind on instruments.

LUBRICATION

Do not use mineral oil or silicone based lubricants as they cannot be easily removed by cleaning. These types of lubricants may coat surface microorganisms and reduce the effectiveness of steam sterilization. They may even provide nutrients for microbial growth.

It is recommended that a water soluble lubricant compatible with the sterilization method be used for lubrication, unless otherwise specified on the labeling.

Follow manufacturer's written instructions for dilution, temperature and whether to rinse or soak the instrument. To protect instruments from staining and rusting during sterilization and storage, they should be lubricated with a water-soluble, preserved lubricant after each cleaning. Since effective ultrasonic cleaning removes all lubricant, relubrication is important. The lubricant should contain a chemical preservative to prevent bacterial growth in the lubricant bath. The bath solution should be made with demineralized water. A lubricant containing a rust inhibitor helps prevent electrolytic corrosion of points and edges. Immediately after cleaning, instrument should be immersed for 30 seconds and allowed to drain off, not wiped off. A lubricant film will remain through the sterilization to protect them during storage.

RESPONSIBILITIES OF THE USER

Health care personnel should wear personal protective equipment like a face mask, goggles and gloves when handling, cleaning or transporting contaminated instruments .

The healthcare facility is responsible for making sure that personnel, equipment and the environment necessary to maintain safe handling and reprocessing of surgical instruments can be achieved. Reference to personnel considerations in ANSI/ AAMI ST79 Comprehensive Guide to Steam Sterilization and Sterility Assurance in Health Care Facilities is recommended.

MANUAL CLEANING INSTRUCTIONS

Instruments must be cleaned and sterilized immediately after surgery. Failure to do so could result in organic soil drying on the instrument thereby making them more difficult to clean and sterilize.

At place of use

Remove all visible soil using a lint free towel. Remove all blood and stains with a non-shedding towel.

Cannulated instruments should be flushed with purified water to remove soil and debris in order to prevent it from drying on the instrument thereby making it challenging to clean. Place instruments in a tray covered with towels dampened with purified water to prevent soil and debris from drying. Clean instruments within 30 minutes of use.

At dedicated cleaning area

NOTE: instruments requiring disassembly should be taken apart prior to cleaning.

Take the covered instruments to the cleaning area.

Rinse instruments under cool water (<104°F/ 40°C) to remove gross soiling (large particles of blood and contamination) within 30 minutes after use.

Do not use a fixating detergent or hot water (>104°F/ 40°C) as this can cause fixation of residual matter.

Prepare an enzymatic solution with tap water.

Follow manufacturer's written instructions for dilution and temperature.

Soak/submerge the instrument in the enzymatic cleaner for at least 45 minutes after a steady stream of bubbles (indicating evacuating air) is no longer observed. Actuate all moving parts during soak while flushing all inside channels with a syringe. Use a soft bristle brush to remove debris, paying special attention to threads, crevices and hard to reach areas. Brush all inside channels with a suitable diameter brush.

Use an appropriately sized brush for lumen being cleaned because if it is too small it will not clean thoroughly and if it is too big it will not fit into the lumen.

Instruments with lumens should be soaked vertically or thoroughly flushed to prevent trapped air bubbles and to ensure complete contact with cleaning solutions.

Rinse and flush instruments under running tap water for 3 minutes.

Use a syringe to flush all hard to reach areas.

Scrub entire instrument including all lumens and cavities.

Scrub until all visible soil (blood and contamination) is removed.

Remove instrument from enzymatic solution and ensure all cavities are flushed in the water stream.

Rinse with running deionized or tap water for 3 minutes.

Prepare a second enzymatic solution with tap water.

Follow manufacturer's written instructions for dilution, temperature and whether to rinse or soak the instrument .

Ultrasonic Cleaner

Submerge instrument in the enzymatic solution and sonicate for at least 45 minutes.

Remove the instrument from the ultrasonic cleaner. Ensure all cavities are flushed. Rinse for 3 minutes.

Rinse with running deionized or tap water for 3 minutes.

Inspect instrument to make sure all blood and contamination have been removed.

Do not proceed with reprocessing a soiled instrument.

Cleaning (Automated): Instrument must be cleaned manually prior to using the automated cleaner. Load instrument so that the instrument can drain.

Drying: Instrument may be dried with a lint free towel. Clean filtered compressed air may be used to dry inside channels and highly inaccessible areas.

Instrument must be properly dried for effective sterilization.

Maintenance of sterility is based on wrap used, wrapping method and storage conditions.

Wrapping method is based on healthcare personnel wrapping technique. Refer to ANSI/AAMI ST79.

Follow wrap manufacturers written instructions for use and sterility maintenance. Wraps used during the steam sterilization process are to be FDA cleared wraps.

The following sterilization instructions have been validated to a Sterility Assurance Level (SAL) of 10^{-6} .

Only legally marketed devices, accessories and chemicals should be used for the reprocessing procedure.

Local laws and ordinances must be followed by the user when they are more stringent than those described in these instructions.

Table 1: Steam Sterilization Cycle Parameters for Medical Facilities In the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	DRYING TIME
Steam	Gravity-Displacement	270°F (132°C)	15 minutes	45 minutes
Steam	Dynamic-Air-Removal	270°F (132°C)	4 minutes	45 minutes

The sterilization cycles listed in Table 2 are not considered by the FDA to be standard sterilization cycles.

To facilitate proper drying open the sterilizer door slightly and leave the devices inside for at least 10 minutes to facilitate proper drying and reduce the chances of forming condensation.

After the drying phase, allow the devices to cool-down for at least 15 minutes on the sterilizer cart outside the steam sterilizer. During this cooling phase the packaged devices should not be touched or transferred to cold metal racks or shelves for cooling as this could cause condensation to form which could contaminate the packaged devices.

Table 2: Steam Sterilization Cycle Parameters for Medical Facilities outside the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	DRYING TIME
Steam	Gravity-Displacement	273°F (134°C)	20 minutes	45 minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	3 minutes	45 minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	4 minutes	45 minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	20 minutes	45 minutes

Instrument must be properly dried for effective sterilization.

NOTE: Chamber size and chamber load differences may exist between industrial and health care facility sterilizer models. The sterilization parameters listed in Tables 1 and 2 can be achieved in both health care facility and larger, industrial sterilizer models.

Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g. temperatures, times) used for their equipment.

EXPLANATION OF SYMBOLS	
 Batch code	 Catalog number
 Authorized European representative	 Product is sold non-sterile
 Only CAUTION: Federal law (USA) restricts these devices to sale by or on the order of a physician.	 This instrument complies with the European Medical Device Directive MDD 93/42/EEC
 Manufacturer	 Consult Instructions for Use

Medical Device

CE mark is only valid if it is also printed on the product and product label.

MDSS GmbH, Schiffgraben 41, 30175 Hannover, Germany

Gauthier Biomedical, Inc., 2221 Washington Street, Grafton, WI 53024 USA

<http://www.gauthierbiomedical.com>

+1-262-546-0010

© 2013-2020 Gauthier Biomedical, Inc. All rights reserved.

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EN

IMPORTANT INFORMATION FOR REUSABLE SURGICAL INSTRUMENTS FROM GAUTHIER BIOMEDICAL.

SCOPE

This instruction for use will apply to reusable manual surgical instruments.

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Rod bending or rod cutting instruments have the specific risks of producing broken bits of material. These instruments could apply very high forces and sometimes bits of material can be produced when rods or other materials are bent or cut. This is why it is hazardous to use them in situ. Any material left in the body of the patient could cause an allergic reaction or an infection.

The physician and rest of healthcare personnel should wear PPE (personal protection equipment) like mask, eye protection (goggles), gloves.

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The health care personnel or patient could become infected from improperly cleaned or sterilized instruments. The patient could also become infected due to surgical complications.

The patient could experience pain, discomfort or abnormal sensations when the instrument is present.

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Remove instrument from packaging then clean and sterilize instrument prior to use in surgery.

Instruments must be cleaned and sterilized prior to returning product to distributor. Health care personnel or the responsible person should include a note to notify the recipient that the instrument has been cleaned and sterilized and therefore is not a biohazard.

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Ensure no parts are missing, loose or damaged. These parts include but are not limited to springs, screws and pins .

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Tips close together properly.

Ratchet teeth are not damaged or missing. If the teeth are not visible, perform a check by connecting a mating instrument and cycle instrument through at least one full rotation in each direction.

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Highly aggressive agents such as Sodium hydroxide (NaOH) or Sodium hypochlorite(NaOCl) commonly known as bleach may discolor or corrode instruments.

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LUBRICANTS

Do not use mineral oil or silicone based lubricants as they cannot be easily removed by cleaning. These types of lubricants may coat surface microorganisms and reduce the effectiveness of steam sterilization. They may even provide nutrients for microbial growth. Oil based lubricants

may be used to lubricate the internal mechanisms of powered instruments if specified in the manufacturer's written instructions.

It is recommended that a water soluble lubricant compatible with the sterilization method be used for lubrication.

Follow manufacturer's written instructions for dilution, temperature and whether to rinse or soak the instrument.

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Place instruments in a tray of water covered with damp towels. Clean instruments within 30 minutes of use.

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Soak/submerge the instrument in an enzymatic cleaner for at least 45-60 minutes after a steady stream of bubbles (indicating evacuating air) is no longer observed. Actuate all moving parts during soak while flushing all inside channels with a syringe. Use a soft bristle brush to remove debris, paying special attention to threads, crevices and hard to reach areas. Brush all inside channels with a suitable diameter brush.

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Use appropriate size brush for lumen being cleaned because if it is too small it will not clean thoroughly and if it too big it will not fit into the lumen.

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Prepare enzymatic solution with tap water.

Follow manufacturer's written instructions for dilution, temperature and whether to rinse or soak the instrument.

Soak/submerge the instrument in an enzymatic cleaner for at least 45-60 minutes after a steady stream of bubbles (indicating evacuating air) is no longer observed. Actuate all moving parts during soak while flushing all inside channels with a syringe. Use a soft bristle brush to remove debris, paying special attention to threads, crevices and hard to reach areas. Brush all inside channels with a suitable diameter brush.

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Brush all inside channels with a suitable diameter brush. Rinse with running deionized or tap water for 3 minutes.

Prepare a second enzymatic solution

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Ultrasonic Cleaner: Submerge instrument in the enzymatic solution and sonicate for 45 - 60 minutes.

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Wraps used during the steam sterilization process are to be FDA cleared wraps.

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Table 1: Steam Sterilization Cycle Parameters for Medical Facilities in the United States and its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	DRYING TIME
Steam	Gravity-Displacement	270°F (132°C)	15 minutes	30 minutes
Steam	Dynamic-Air-Removal	270°F (132°C)	4 minutes	30 minutes

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Table 2: Steam Sterilization Cycle Parameters for Medical Facilities outside the United States and its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	DRYING TIME
Steam	Gravity-Displacement	273°F (134°C)	20 minutes	30 minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	4 minutes	30 minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	20 minutes	30 minutes

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<http://www.gauthierbiomedical.com>

+1-262-546-0010

© 2013 Gauthier Biomedical, Inc. All rights reserved.

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





**Tecomet
Kenosha, WI 53144 USA
Toll Free 855-925-4500
262-925-4500**

KLIFU0002-06

Reprocessing Surgical Instruments

Issue Date: 2021/11/10

Product

Surgical instruments

Warning Notice

Instruments are provided non-sterile. Clean and sterilize prior to each use. Disassemble instruments with removal parts prior to cleaning and sterilization. During cleaning of tight lumina or drill holes special attention might be required. **USA – RX ONLY – Caution: Federal (USA) law restricts this device to sale by or on the order of a licensed healthcare practitioner.**

Advise

Devices which cannot be reprocessed by the following processes are labeled accordingly. Individual procedures accompanying those devices must be followed. In case of damage, the device should be reprocessed before returning to the manufacturer for repair.

Reprocessing Limitations:

Repeated reprocessing has minimal effect upon manual instruments. End of life is determined by wear and damage due to use.

Preparation at the Point of Use:

Remove gross soiling by submerging the instrument into cold water (<40°C) immediately after use.

Don't use a fixating detergent or hot water (>40°C) as this can cause the fixation of residua which may influence the result of the reprocessing process. Remove surface contamination with paper tissue.

Containment / Transportation:

Follow hospital protocols when handling contaminated and

bio-hazardous materials. Instrument should be cleaned within 30 minutes after use to minimize the potential of staining, damage, and drying.

Preparation for Decontamination:

If possible, the devices must be reprocessed in a disassembled or opened state.

Pre-Cleaning:

Not required.

Cleaning:

Manual Cleaning Process:

- Rinse under running tap water (<40°C) until all visible soil has been removed. If needed, a soft bristle brush should be used to remove visible soil.
- Submerge instruments in an enzymatic detergent (If ultrasonic bath is used, ultrasonic process of 3 minutes and ultrasonic frequency of 35 kHz have been shown to be effective). Follow detergent manufacturer's instructions.
- Rinse the instrument under running tap water to remove the detergent.

Drying:

Manual Drying:

- Dry the instrument with a lint-free towel. The instrument may never be heated up >140°C. To avoid water residues, insufflate cavities of instruments by using sterile compressed air.

Inspection / Functional Testing:

1. Carefully inspect each device to ensure that all visible blood and soil has been removed.
2. Visually inspect for damage and / or wear.

3. Check the action of moving parts (such as hinges and box-locks) to ensure smooth operation throughout the intended range of motion.
4. Check instruments with long slender features (particularly rotating instruments) for distortion.
5. Where instruments form part of a larger assembly, check that the devices assemble readily with mating components.

Note: If damage or wear is noted that may compromise the function of the instrument, do not use the instrument, and notify the appropriate responsible person.

Packaging:

Appropriate packaging for sterilization.

Sterilization:

Sterilize instruments by applying a fractionated pre-vacuum process under consideration of the respective country requirements.

We have validated and recommend using the following parameters for the fractionated pre-vacuum process or gravity-process: 132°C (270°F) for 10 minutes and 20 minutes drying time, 121°C (250°F) for 30 minutes and 15 to 30 minutes drying time.

Flash sterilization is not allowed on lumen instruments!

Warning:

Tecomet recommends that these devices are not compatible with automated cleaning processes.

Storage:

Storage of sterilized instruments in a dry, clean, and dust-free environment at modest temperatures of 5°C to 40°C.

Maintenance:

Lubricate hinges, threads and other moving parts with a commercial water-based surgical grade instrument lubricant (such as instrument milk) to reduce friction and wear. Follow lubricant manufacturer's instructions.

Additional Advices:

It is the duty of the user to ensure that the reprocessing processes are followed; that resources and materials are available to capable personnel; and that hospital protocols are followed. State of the art and often national law requiring these processes and included resources to be validated and maintained properly.

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2021-11-10

Medtronic
ENGLISH
AUSTRALIAN SPONSOR:

Medtronic Australasia Pty Ltd
 97 Waterloo Rd
 North Ryde, NSW 2113
 Australia

Medtronic Sofamor Danek USA, Inc.
 1800 Pyramid Place
 Memphis, TN 38132
 Telephone 800 933 2635 (In U.S.A.)
 901 396 3133 (Outside U.S.A.)
 Fax 901 396 0356

ECIREP

Medtronic B.V.
 Earl Bakkenstraat 10
 6422 PJ Heerlen
 The Netherlands
 Tel: + 31 45 566 80 00

IMPORTANT INFORMATION ON THE POWEREASE™ SYSTEM WORKING ENDS
PURPOSE

The instruments are working ends that are intended to be used in surgical procedures in the thoracic, lumbar, and sacral spine.

DESCRIPTION

The working ends are compatible with the IPC® POWEREASE™ System and the NIM-ECLIPSE™ System. The working ends consist of instruments such as taps, drill bits, reduction nut drivers, screw drivers, rod cutter, post cutter, set screw break-off tool, and sleeves. The working ends have a manual alternative.

The working ends are compatible with the CD HORIZON® SOLERA® Spinal System and the TSRH® 3Dx Spinal System.

The taps, screw drivers, drill bits, and sleeves are also compatible with Medtronic's NIM-ECLIPSE™ Spinal System.

INDICATIONS

The IPC® POWEREASE™ System is indicated for drilling, tapping, and driving screws and working end attachments during spinal surgery, including open and minimally invasive procedures. The IPC POWEREASE™ System is also used in the placement or cutting of screws, posts, and rods.

DO NOT IMPLANT THE INSTRUMENTS.

If there is any doubt or uncertainty concerning the proper use of this instrument, please contact MEDTRONIC. Any available surgical techniques will be provided at no charge.

MEDTRONIC does not and cannot warrant the use of this instrument nor any of the component parts upon which repairs have been made or attempted, except as performed by MEDTRONIC or an authorized MEDTRONIC repair representative. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded. See the MDT Catalog or price list for further information about warranties and limitations of liability.

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See the IPC POWEREASE™ System package insert for instructions to assemble the working ends to the POWEREASE™ driver.

Taps, Drill Bits and Drivers

Once a tap, reduction nut driver or drill bit is attached to the POWEREASE™ Driver, the instrument is ready to use.


- To attach an MAS, FAS, or RMAS screw to the screwdriver working end, align the driving tip within the screw and tighten the screwdriver's sleeve into the screw.
- To attach a Post Screw to the Post Screwdriver, first assemble the screwdriver base components by threading them together.
 - Insert the screw into the tip so that the screwdriver captures the screw rotationally.
- To attach a Threaded Post Screw to the Threaded Post Screwdriver, insert the screw into the tip so that the screwdriver captures the screw rotationally.
- To attach an MPA screw to the Locking MPA™ Screwdriver, first assemble the screwdriver base components by sliding the Locking MPA™ Driver Shaft into the Locking MPA™ Sleeve and threading them together.
 - Insert the screw into the tip so that the screw is fully seated.
 - Rotate the cap clockwise on the Locking MPA™ Sleeve until the Locking MPA™ Screwdriver securely holds the MPA™ Screw.
 - To lock the cap, pull down and rotate the cage on the Locking MPA™ Sleeve, then release so that the teeth between the cage and the cap engage.

Instructions for Manual Use

- Insert working end into the ratcheting handle provided in the implant set by aligning the working end with the ratcheting handle.
- Insert until fully seated within the handle.

Note

- Only use ratcheting handles compatible with the specific working end.
- Drill bits are single use only.
- Taps and drill bits are marked in 1cm increments. The depth indication is measured from the tip of the working end to the horizontal band.

Optional

A sleeve is provided which can be used as a secondary handle for stability while tapping the bone or driving a screw with the POWEREASE™ Driver.

1. Determine the proper size sleeve to use with the tap, drill bit, or screwdriver.
2. Slide the sleeve onto the working end so that the wide part of the sleeve is at the back of the working end. The sleeve will hit a physical stop. For taps 8mm and larger and MPA™ Screwdriver, the 8.0mm sleeve must slide over the back of the working end prior to attachment to the POWEREASE™ Driver.
3. Rotate the colored Cam on the sleeve in the direction of the arrow.

Set Screw Break Off (SSBO) Instrument

1. Select the appropriate Counter Torque.
2. Attach the Counter Torque to the SSBO Instrument by sliding the Counter Torque over the SSBO Breaker Shaft.
3. Press the lock button and hold it.
4. Slide the Counter Torque into the Body of the SSBO Instrument.
5. Release the Lock Button.
6. Attach the SSBO Instrument to the POWEREASE™ Driver.
7. Align the SSBO Instrument Breaker Shaft with the set screw and press down to seat the Breaker Shaft fully over the set screw.
8. Once the SSBO Instrument Breaker Shaft is fully seated onto the set screw, press the SSBO instrument down and rotate the SSBO Instrument until the Counter Torque cut-outs are aligned with the rod.
9. Once the Counter Torque is aligned with the rod, fully seat the Counter Torque over the rod by pressing down on the SSBO Instrument.
10. Once the Counter Torque is fully seated, run the POWEREASE™ Driver in forward mode until the set screw is broken. The broken set screw heads are stored in the Breaker Shaft.
11. Pull the SSBO instrument away from the construct.

Instructions for Purging Stored Set Screw Heads from the Breaker Shaft

The Breaker Shaft of the SSBO Instrument can hold between 10 and 14 broken set screw heads depending on the spinal system used. When the Breaker Shaft becomes full of broken set screws, or at the completion of surgery, it is necessary to purge the broken set screw heads.

1. Remove the Counter Torque from the SSBO Instrument.
2. Insert an Obturator into the front of the Breaker Shaft of the SSBO instrument.
3. Press inward until all the broken set screw heads are ejected from the side port of the Breaker Shaft.

Note

- Light rotation of the SSBO Instrument may be necessary to seat onto the set screw.
- If the Counter Torque resists being fully seated, briefly press the POWEREASE™ Driver trigger.
- The SSBO Instrument's drive train is not connected to the Breaker Shaft until the Counter Torque is fully seated onto the rod.

Instructions for Manual Use

1. Follow all the steps described for attaching the SSBO Instrument to the POWEREASE™ Driver, with the exception of Step 2.
2. For Step 2, instead of the POWEREASE™ Driver, a manual handle with the Square Adapter attached may be used to drive the SSBO Instrument.

Post Cutter

The Post cutter is used to cut long post screws that are made from Titanium Alloy and have a diameter of 5.2mm. If using a break-off set screw, the set screw head must be removed prior to using the Post Cutter.

1. Attach the Post Cutter to the POWEREASE™ Driver.
2. Ensure that the Post Cutter is in the open position by running the Post Cutter in reverse until the Post Cutter mechanism is no longer changing, as indicated by the horizontal lines at the pivot point on the Post Cutter arms.
3. Position the Post Cutter over the post so the Post Cutter is flush with the connector.
4. Run the Post Cutter in forward until the post is cut.

To Remove the Cut Portion of the Post:

5. Move the Post Cutter away from the patient.
6. Run the Post Cutter in reverse. The cut post will be able to slide out the back of the Post Cutter.

Instructions for Manual Use

In step 1 of the previous section, attach the Post Cutter to a manual handle with the Square Adapter instead of the POWEREASE™ Driver. Follow the remaining steps for using the Post Cutter.

Rod Cutter

The Rod Cutter is used to cut surgical rods of material types such as Titanium Alloy, Titanium CP, CHROMALLOY™, and CHROMALLOY™ Plus in diameters ranging from 4.75mm to 6.0mm.

1. Attach the Rod Cutter to the POWEREASE™ Driver.
2. Ensure that the Rod Cutter is in the open position by running the Rod Cutter in reverse until the Rod Cutter mechanism is no longer changing, as indicated by the horizontal through holes being aligned.
3. Ensure that both ends of the rod are tightly held to reduce the risk caused by projectile rod pieces.

NOTE: Do not cut rods in-situ.

4. There are two methods of cutting a rod: Scissor Cutting and Shear Cutting.
 - a. For Scissor Cutting, place the rod between the cutting tips so that the rod is touching both surfaces.
 - b. For Shear Cutting, pass the rod through the hole that correlates to its size (5.5mm rod through a 6.0mm hole). The cutting point is 6mm beyond the entry holes.
5. Run the Rod Cutter in forward until the rod is cut. It may be necessary to run the Rod Cutter in reverse to remove the cut portion of the rod.

Instructions for Manual Use

Follow all the steps described in the previous section for using the Rod Cutter with the POWEREASE™ Driver, with the exception of step 1. In the place of the POWEREASE™ Driver, a manual handle with the Square Adapter attached may be used to drive the Rod Cutter.

Note for Post Cutter and Rod Cutter

The Post Cutter and the Rod Cutter should be lubricated between every surgical use to ensure smooth performance and a long service life. Lubricate hinges, threads, and other moving parts with a commercial water-based surgical grade instrument lubricant (such as instrument milk) to reduce friction and wear. Follow lubricant manufacturer's instructions.

WARNINGS

- Breakage, slippage, misuse, or mishandling of instruments, such as on sharp edges, may cause injury to the patient or operative personnel.
- Improper maintenance, handling, or poor cleaning procedures can render the instrument unsuitable for its intended purpose or even dangerous to the patient or surgical staff.
- There are particular risks involved in the use of instruments used for bending and cutting rods. The use of these types of instruments can cause injury to the patient by virtue of the extremely high forces which are involved. **DO NOT CUT RODS IN SITU.** In addition, any breakage of an instrument or the implant in this situation could be extremely hazardous.
- It is important that the surgeon exercise extreme caution when working in close proximity to vital organs, nerves, or vessels and that the forces applied while correcting the position of the instrumentation is not excessive, such that it might cause injury to the patient.

PRECAUTIONS

- Excessive force applied by instruments to implants can dislodge devices, particularly hooks.
- Never expose instruments to temperatures in excess of 135 °C that may considerably modify the physical characteristics.
- Extreme care should be taken to ensure that this instrument remains in good working order. During the procedure, successful utilization of this instrument is extremely important. Instruments should not be bent or damaged in any way. Misuse of instruments resulting in corrosion, "freezing-up", scratching, loosening, bending, or fracture of any or all sections of an instrument may inhibit or prevent proper function.
- These instruments should be carefully placed on trays, cleaned after each use, and stored in a dry environment.
- Do not use this instrument for any action for which it was not intended.
- Regularly review the operational state of all instruments and, if necessary, make use of repair and replacement services.
- To avoid injury, the instrument should be carefully examined for functionality or damage prior to use. A damaged instrument should not be used. Additional back-up instruments should be available.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results.
- Proper patient selection and operative care are critical to the success of the surgery and avoidance of injury during surgery. Read and follow all other product information supplied by the manufacturer of the implants or the instruments.
- Special precautions are needed during pediatric use. Care should be taken when using instruments in pediatric patients since these patients can be more susceptible to the stresses involved in their use.
- Some surgeries require the use of instruments which incorporate a measuring function. Ensure that these are not worn and any surface engravings are clearly visible.

POSSIBLE ADVERSE EFFECTS

- Nerve damage, paralysis, pain, or damage to soft tissue, visceral organs, or joints.
- Infection if instruments are not properly cleaned and sterilized.
- Pain, discomfort, or abnormal sensations resulting from the presence of the instrument.
- Nerve damage due to surgical trauma.
- Dural leak in cases of excessive load application.
- Impingement of close vessels, nerves, and organs by slippage or misplacement of the instrument.
- Damage due to spontaneous release of clamping devices or spring mechanisms of certain instruments.
- Cutting of skin or gloves of operating staff.
- Bony fracture in cases of deformed spine or weak bone.
- Tissue damage to the patient, physical injury to operating staff, and/or increased operating time that may result from the disassembly of multi-component instruments occurring during surgery.
- The methods of use of instruments are to be determined by the user's experience and training in surgical procedures. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results.

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Instruments supplied sterile are designed for single patient use only. Do not reuse, reprocess, or re-sterilize these products. Reuse, reprocessing, or re-sterilization may compromise the structural integrity of the instrument and/or create a risk of contamination of the instrument, which could result in patient injury, illness, or death.

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PACKAGING

Packages for components should be intact upon receipt. All sets should be checked for completeness and all components should be checked for signs of damage prior to use. Damaged packages or products should not be used and should be returned to MEDTRONIC.

Remove all packaging material prior to sterilization. Only sterile instruments should be used in surgery. Always immediately clean and re-sterilize all instruments used in surgery. Instruments should be thoroughly cleaned prior to re-sterilization. This process must be performed before handling, or before returning product to MEDTRONIC.

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Examination should be thorough and must include a visual and functional inspection of the working surfaces, pivots, racks, spring or torsional operation, cleanliness of location holes or cannulations, and the presence of any cracks, bending, deformation, or distortion, and that all components are complete.

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Visual Inspection

Make certain of the following:

- Laser etchings, engravings, and other markings are legible.
- No cracks are present in instrument handles or any part of the instrument.
- Discoloration, corrosion, stains, or rust do not exist. If present, attempt to remove in accordance with the cleaning instructions provided in the Reprocessing section of this document.
- There is no handle/shaft separation, and the handle-to-shaft connection is secure.
- No cuts or gouges in silicone are present.
- There is no damage (cuts, tears, etc.) to insulation.
- There is no damage to the working ends or tips. The working ends should be free of cracks, sharp edged gouges, and other damage. When applicable, the working ends should be sharp.
- There is no damage to threads.
- All parts are present and free of damage and deterioration. Examples of parts that may be missing, loose, or damaged include set screws, springs, curved springs, pins, and prongs.
- Mating ends are free of damage (nicks, gouges, bends, etc.) that would interfere with the mating function.
- Cannulated instruments with a guide wire or other insertion tool are visually checked.

Functional Inspection

Make certain of the following:

- The parts intended to move will do so freely without sticking, binding, or grinding.
- Springs return the handle of the instrument to its original position.
- Retention tabs hold appropriate mating parts and are not damaged.
- The instrument will function as intended with the appropriate mating parts.
- Ball detents will hold mating parts and are free from damage.
- Sharp edges are sharp to the touch and are not dull, have no nicks, or any other damage.
- Tips meet when appropriate.
- Ratcheting mechanisms are functional. This includes handles, latches, and other mechanisms. All teeth should be present and functional.
- Driver tips are not worn beyond functional use. If necessary, mate the instrument with the appropriate part.

REPROCESSING

Refer to the Reprocessing Instructions for the POWEREASE™ System Working Ends package insert (M7083488184) for the cleaning instructions and sterilization parameters and requirements.

FURTHER INFORMATION

In case of complaint, or for supplementary information, contact MEDTRONIC.

MRI INFORMATION

These instruments should not be introduced into the MR environment.

PRODUCT COMPLAINTS

Any health care professional (e.g., user of this system of products) having any complaints or having experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, or performance should notify the distributor or MEDTRONIC. Further, if any of the implanted spinal system component(s) ever "malfunctions" (i.e., does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any MEDTRONIC product ever "malfunctions" and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax, or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint, and notification of whether or not a written report from the distributor is requested.

© 2013 MEDTRONIC SOFAMOR DANEK USA, Inc. All rights reserved.

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Reprocessing Instructions for the POWEEASE™ System Working Ends

M708348B184 Rev. A

REPROCESSING INSTRUCTIONS FOR THE POWEEASE™ SYSTEM WORKING ENDS
INSTRUCTIONS DE RETRAITEMENT DES EXTRÉMITÉS DE TRAVAIL DU SYSTÈME POWEEASE™
INSTRUCCIONES DE REPROCESAMIENTO PARA LOS EXTREMOS FUNCIONALES DEL SISTEMA POWEEASE™
ANLEITUNGEN ZUR WIEDERAUFBEREITUNG DES POWEEASE™-SYSTEMS ARBEITSSENDEN
ISTRUZIONI SUL RITRATTAMENTO DEL SISTEMA POWEEASE™ ESTREMITÀ DI LAVORO
HERVERWERKINGSINSTRUCTIES VOOR HET POWEEASE™ SYSTEEM WERKUIEINDEN
GENBEARBEJDNINGSANVISNINGER TIL POWEEASE™ SYSTEMET ARBEJDSENDER
BEARBETNINGSINSTRUKTIONER FÖR POWEEASE™-SYSTEMETS ARBETSÄNDAR
INSTRUKSJONER FOR GJENBRUKSBEHANDLING AV POWEEASE™-SYSTEMET ARBEJDSENDER
JÄLLEENKÄSITTELYOHJEET POWEEASE™-JÄRJESTELMÄÄ VARTEN TYÖPÄÄT
ΟΔΗΓΙΕΣ ΕΠΑΝΕΠΕΞΕΡΓΑΣΙΑΣ ΓΙΑ ΤΟ ΣΥΣΤΗΜΑ POWEEASE™ SYSTEM ΛΕΙΤΟΥΡΓΙΚΑ ΑΚΡΑ
INSTRUÇÕES DE REPROCESSAMENTO DO SISTEMA POWEEASE™ EXTREMIDADES DE TRABALHO
POKYNÝ PRO OPAKOVÁNÉ ZPRACOVÁNÍ PRO SYSTÉM POWEEASE™ PRACOVNÍ KONCE
SÜSTEEMI POWEEASE™ TÖÖOTSIKUTE ÜMBERTÖÖTAMISJUHISED
ÚJRAFELDOLGOZÁSI UTASÍTÁSOK A POWEEASE™ RENDSZERHEZ MUNKAVÉGEK
ATTÍRÍŠANAS NORÁDÍJUMI POWEEASE™ SISTÉMAI DARBA UZGAĻI
NURODYMAI DĖL „POWEEASE™“ SISTEMOS APDOROJIMO DARBINIAI ANT GALIAI
INSTRUKCJE PRZYGOTOWANIA DO PONOWNEGO UŻYCIA KOŃCÓWEK ROBOCZYCH SYSTEMU POWEEASE™
POKYNÝ PRE OPÄTOVNÉ SPRACOVANIE SYSTÉMU POWEEASE™ PRACOVNÉ KONCE

2012-01

Medtronic

AUSTRALIAN SPONSOR:

Medtronic Australasia Pty Ltd
97 Waterloo Rd
North Ryde, NSW 2113
Australia

Medtronic Sofamor Danek USA, Inc.
1800 Pyramid Place
Memphis, TN 38132
Telephone 800 933 2635 (In U.S.A.)
901 396 3133 (Outside U.S.A.)
Fax 901 396 0356

Medtronic B.V.
Earl Bakkenstraat 10
6422 PJ Heerlen
The Netherlands
Tel: +31 45 566 8000

ENGLISH

REPROCESSING INSTRUCTIONS FOR THE POWEEASE™ SYSTEM WORKING ENDS

GENERAL CONSIDERATIONS

It is the purpose of this document to provide cleaning and sterilization procedures for the POWEEASE™ System Working Ends.

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Thorough Cleaning of Instruments

It is critical that reusable instruments are thoroughly cleaned after each use following the specific cleaning steps listed in this document. Ineffective or incomplete cleaning can hinder subsequent sterilization activities.

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Cleaning and Rinsing Water

Tap water should be used. De-ionized water should be used for the final rinse step to prevent mineral deposits on instrument surfaces.

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MANUAL CLEANING INSTRUCTIONS

Proper Handling After Use

DO NOT allow the instruments to dry after use and prior to cleaning. Cleaning and subsequent sterilization may be hindered when blood or bloody solutions are allowed to dry on instruments.

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1. Immediately transport the tray containing covered instruments to a work area dedicated to further reprocessing.
2. Rinse instruments under running tap water for 3 minutes.
3. Scrub instruments with an appropriate soft bristle brush to aid in the removal of visible soil.
4. Repeat steps 2 and 3 until all visible soil is removed.
5. Using tap water, prepare an enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures.
6. Place instruments in the enzymatic cleaner, completely submerged, and soak for 45-60 minutes.
7. Remove instruments from the enzymatic cleaner and rinse under running tap water for 3 minutes.
8. Using tap water, prepare a second enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures in an appropriately-sized sonicator.
9. Place instruments in the enzymatic cleaner, completely submerged, and sonicate for 45-60 minutes.
10. Remove instruments from the sonicator and rinse using running tap water for 3 minutes.
11. Repeat rinsing step 10, using de-ionized water for an additional 3 minutes.
12. Dry the instruments using clean, absorbent, non-shedding wipes.
13. Carefully inspect instruments, including any lumens and cavities, to ensure that all contamination has been removed. If any soil is still present, repeat the cleaning process or contact your Medtronic distributor immediately to arrange for disposal or replacement. Do not proceed with reprocessing of a soiled instrument.

Inspection

1. Inspect all instruments prior to use. Carefully check for signs of corrosion, discoloring, pitting, or any other signs of wear.
2. Ensure that the instruments do not appear bent or damaged. Also inspect the handles for discoloration, cuts, or other signs of damage.

SANITIZATION

The Working Ends must be sterilized prior to initial use, or in adherence to these reprocessing instructions before re-use. Sanitization with disinfectant solutions or chemicals is unnecessary and not recommended.

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Ethylene oxide (EO), gas plasma, gamma irradiation, chemical vapor, or dry heat sterilization methods are NOT recommended for sterilization. Steam/moist heat is the recommended method of sterilization.

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Sterilization wraps, pouches, chemical indicators, biological indicators, and sterilization trays should be FDA-cleared for the selected sterilization cycle parameters.

The sterilization instructions below have been validated to a sterility assurance level of 10^{-4} .

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The sterilization cycles listed in Table 2 are not considered by the Food and Drug Administration (FDA) to be standard sterilization cycles. It is the user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

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PRODUCT COMPLAINTS

Any health care professional (e.g., user of this system of products) having any complaints or having experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, or performance, should notify the distributor or MEDTRONIC. Further, if any of the implanted spinal system component(s) ever "malfunctions" (i.e., does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any MEDTRONIC product ever "malfunctions" and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax, or written correspondence. When filing a complaint, please provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint, and notification of whether or not a written report from the distributor is requested.

©2012 MEDTRONIC SOFAMOR DANEK USA, Inc. All rights reserved.

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Medtronic NON-REUSABLE and NON-STERILE INSTRUMENTS

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
ENGLISH

IMPORTANT INFORMATION FOR MEDTRONIC NON-REUSABLE AND NON-STERILE INSTRUMENTS

DESCRIPTION

Non-reusable and non-sterile instruments are made out of a variety of materials commonly used in orthopedic and neurological procedures which meet available national or international standards specifications.

None of these instruments should ever be reused under any circumstances.

INTENDED USE

These orthopedic manual surgical instruments are intended for use in surgical procedures to manipulate tissue, bone, or for use with other devices in orthopedic surgery. An instrument may incorporate a measuring function, which has uses as described on the label and the instrument.

DO NOT IMPLANT THE INSTRUMENTS.

Medtronic does not and cannot warrant the use of these instruments nor any of the component parts upon which repairs have been made or attempted, except as performed by Medtronic or an authorized Medtronic repair representative. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

WARNINGS

- These instruments should never be reused, reprocessed, or resterilized under any circumstances. Reuse, reprocessing, or resterilization may compromise the structural integrity of these instruments and create a risk of contamination, which could result in patient injury, illness, or death.
- Breakage, slippage, misuse, or mishandling of instruments, such as on sharp edges, may cause injury to the patient or operative personnel.
- Improper maintenance, handling, or poor cleaning procedures can render instruments unsuitable for intended purposes or even dangerous to the patient or surgical staff.
- There are particular risks involved in the use of instruments used for bending and cutting rods. The use of these types of instruments can cause injury to the patient by virtue of the extremely high forces which are involved. Do not cut rods in situ. In addition, any breakage of an instrument or the implant in this situation could be extremely hazardous. The physical characteristics required for many instruments do not permit them to be manufactured from implantable materials, and if any broken fragments of instruments remain in the body of a patient, the patient could have allergic or infectious consequences.
- It is important that the surgeon exercise extreme caution when working in close proximity to vital organs, nerves or vessels, and that the forces applied while correcting the position of the instrumentation is not excessive, such that it might cause injury to the patient.

PRECAUTIONS

- Excessive force applied by instruments to implants can dislodge devices, particularly hooks.
- Never expose instruments to temperatures in excess of 135°C. This may considerably modify the physical characteristics.
- Extreme care should be taken to ensure instruments are in good working order. During the procedure, successful utilization of instruments is extremely important.
- Instruments should not be bent or damaged in any way. Misuse of instruments, resulting in corrosion, "freezing-up", scratching, loosening, bending, or fracture of any or all sections of an instrument may inhibit or prevent proper function.
- Instruments should be carefully placed on trays and discarded after use.
- Do not use instruments for any action for which they are not intended.
- To avoid injury, instruments should be carefully examined for functionality or damage prior to use. Damaged instruments should not be used. Additional back-up instruments should be available.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results.
- Proper patient selection and operative care are critical to the success of the surgery and avoidance of injury during surgery.
- Read and follow all other product information supplied by the manufacturer of the implants or the instruments.
- Special precautions are needed during pediatric use. Care should be taken when using instruments in pediatric patients, since these patients can be more susceptible to the stresses involved in their use.
- Some surgeries require the use of instruments that incorporate a measuring function. Ensure that these are not worn and any surface engravings are clearly visible.

M708348B359-B

POSSIBLE ADVERSE EFFECTS

- Nerve damage, paralysis, pain, or damage to soft tissue, visceral organs, or joints.
- Infection if instruments are not properly cleaned and sterilized.
- Pain, discomfort, or abnormal sensations resulting from the presence of the device.
- Nerve damage due to surgical trauma.
- Dural leak in cases of excessive load application.
- Impingement of close vessels, nerves, and organs by slippage or misplacement of instruments.
- Damage due to spontaneous release of clamping devices or spring mechanisms of certain instruments.
- Cutting of skin or gloves of operating staff.
- Bony fracture in cases of deformed spine or weak bone.
- Tissue damage to the patient, physical injury to operating staff, and/or increased operating time that may result from the disassembly of multi-component instruments occurring during surgery.
- The methods of use of instruments are to be determined by the user's experience and training in surgical procedures. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results.

PHYSICIAN NOTE: Although the physician is the learned intermediary between the company and patient, the important medical information given in this document should be conveyed to the patient.

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Remove all packaging material prior to sterilization. Only instruments that have been sterilized should be used in surgery. Always immediately discard all instruments used in surgery.

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Never use instruments with obvious signs of excessive wear, damage, or that are incomplete or otherwise non-functional.

Visual Inspection

Make certain of the following:

- Laser etchings, engravings, and other markings are legible.
- No cracks are present in instrument handles or any part of the instrument.
- Discoloration, corrosion, stains, or rust does not exist.
- There is no handle/shaft separation, and the handle-to-shaft connection is secure.
- No cuts or gouges in silicone are present.
- There is no damage (cuts, tears, etc.) to the insulation.
- There is no damage to the working ends or tips. The working end should be free of cracks, sharp edged gouges, and other damage. When applicable, the working end should be sharp.
- There is no damage to threads.
- All parts are present and free of damage. Examples of parts that may be missing, loose, or damaged include set screws, springs, curved springs, pins, and prongs.
- Mating ends are free of damage (nicks, gouges, bends, etc.) that would interfere with the mating function.
- Cannulated instruments with a guide wire or other insertion tool are visually checked.

Functional Inspection

Make certain of the following:

- The parts intended to move will do so freely, without sticking, binding, or grinding.
- Springs return the handle of the instrument to its original position.
- Retention tabs hold appropriate mating parts and are not damaged.
- The instrument will function as intended with the appropriate mating parts.
- Ball detents will hold mating parts and are free from damage.
- Sharp edges are sharp to the touch and are not dull, have no nicks, or any other damage.
- Tips meet when appropriate.
- Ratcheting mechanisms are functional. This includes handles, latches, and other mechanisms. All teeth should be present and functional.
- Driver tips are not worn beyond functional use. If necessary, mate instruments with appropriate parts.

SANITIZATION

Medtronic non-reusable and non-sterile instruments are considered critical devices and must be sterilized prior to initial use. Sanitization with disinfectant solutions or chemicals is unnecessary and not recommended.

STERILIZATION INSTRUCTIONS

Instruments must be sterilized prior to initial use. Ethylene oxide (EO), gas plasma, gamma irradiation, chemical vapor, or dry heat sterilization methods are NOT recommended for sterilization of these instruments. Steam/moist heat is the recommended method of sterilization.

The following sterilization instructions have been validated to a sterility assurance level of 10⁻⁶.

1. Wrap or pouch the parts.
2. Inspect the packaging to ensure no rips, punctures, or seal failures are present in or on the packaging prior to loading into the sterilizer.
3. Load the parts into the sterilizer by following the sterilizer manufacturer's recommended loading procedures and load configurations.
4. Follow the sterilizer manufacturer's recommended procedures to program the sterilizer with any one of the sets of sterilization cycle parameters found in Table 1 (for medical facilities in the US and its territories) and Table 2 (for medical facilities outside the US and its territories).

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It is the user's responsibility to use only sterilizers and accessories (i.e., sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) cleared by the Food and Drug Administration (FDA) for the selected sterilization cycle specifications (time and temperature). The sterilization cycles listed in Table 2 are not considered by the FDA to be standard sterilization cycles.

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Table 2: Sterilization Cycle Parameters for Medical Facilities outside the United States and Its Territories

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INSPECTION AND REASSEMBLY

1. Carefully inspect the instrument for damage by performing a thorough examination as referred to the Examination section of this document.
2. If the instrument is damaged, contact your Medtronic distributor immediately to arrange for disposal or replacement. Do not proceed with a damaged device.
3. If the instrument was disassembled prior to sterilization, reassemble.

PRODUCT COMPLAINTS

Any health care professional (e.g., customer or user of this system of products) who has any complaints or who has experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, and/or performance should notify the distributor (Medtronic). Further, if any of the implanted spinal system component(s) ever "malfunctions" (i.e., does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any Medtronic product ever malfunctions and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, FAX, or written correspondence. When filing a complaint, provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint, and notification of whether or not a written report from the distributor is requested.

MRI INFORMATION

These instruments have not been evaluated for safety, heating, migration, or compatibility in the magnetic resonance (MR) environment.

FURTHER INFORMATION

If there is any doubt or uncertainty concerning the proper use of this instrument, contact Medtronic. Any available surgical techniques will be provided at no charge.

© 2013 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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ENGLISH

IMPORTANT INFORMATION ON THE CD HORIZON™ SPINAL SYSTEM

PURPOSE

The CD HORIZON™ Spinal System is intended to help provide immobilization and stabilization of spinal segments as an adjunct to fusion of the thoracic, lumbar, or sacral spine.

DESCRIPTION

The CD HORIZON™ Spinal System consists of a variety of shapes and sizes of rods, hooks, screws, CROSSLINK™ Plates, staples, and connecting components, as well as implant components from other Medtronic spinal systems which can be rigidly locked into a variety of configurations, with each construct being tailor-made for the individual case.

A subset of CD HORIZON™ Spinal System components may be used for posterior pedicle screw fixation in pediatric cases. These constructs may be comprised of a variety of shapes and sizes of rods (ranging in diameter from 3.5mm to 6.35mm), hooks, screws, CROSSLINK™ Plates and connecting components. Similar to the CD HORIZON™ implants used in adult cases, these components can be rigidly locked into a variety of configurations, with each construct being tailor-made for the individual case.

Certain components within the CD HORIZON™ Spinal System are specifically excluded for use in pediatric patients. These include PEEK rods, Shape Memory Alloy Staples, SPIRE™ Plates, and DYNALOK™ bolts. All screws used in pediatric cases are only cleared for use via a posterior approach. All of the components used in pediatric cases are fabricated from medical grade stainless steel, medical grade titanium, titanium alloy, and medical grade cobalt-chromium-molybdenum alloy.

Certain implant components from other Medtronic spinal systems can be used with the CD HORIZON™ Spinal System in non-pediatric cases. These components include TSRH™ rods, hooks, screws, plates, CROSSLINK™ plates, connectors, staples and washers, GDLH™ rods, hooks, connectors and CROSSLINK™ bar and connectors; LIBERTY™ rods and screws; DYNALOK™ PLUS and DYNALOK CLASSIC™ bolts along with rod/bolt connectors; and Medtronic Multi-Axial rods and screws. Note that certain components are specifically designed to connect to specific rod diameters, while other components can connect to multiple rod diameters. Care should be taken so the correct components are used in the spinal construct.

CD HORIZON™ hooks are intended for posterior use only. CD HORIZON™ staples and CD HORIZON™ ECLIPSE™ rods and associated screws are intended for anterior use only. However, for patients of smaller stature and pediatric patients, CD HORIZON™ 4.5mm rods and associated components may be used posteriorly.

The CD HORIZON™ Spinal System implant components are fabricated from medical grade stainless steel, medical grade titanium, titanium alloy, medical grade cobalt-chromium-molybdenum alloy, or medical grade PEEK OPTIMA-LT1. Certain CD HORIZON™ Spinal System components may be coated with hydroxyapatite. No warranties, expressed or implied, are made. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

Never use stainless steel and titanium implant components in the same construct.

Medical grade titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy may be used together. Never use titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy with stainless steel in the same construct.

The CD HORIZON™ Spinal System also includes anterior staples made of Shape Memory Alloy (Nitinol – NiTi). Shape Memory Alloy is compatible with titanium, titanium alloy, and cobalt-chromium-molybdenum alloy. Do not use with stainless steel. These staples are not to be used in pediatric patients.

PEEK OPTIMA-LT1 implants may be used with stainless steel, titanium or cobalt-chromium-molybdenum alloy implants. CD HORIZON™ PEEK Rods are not to be used with CROSSLINK™ Plates or in pediatric patients.

To achieve best results, do not use any of the CD HORIZON™ Spinal System implant components with components from any other system or manufacturer unless specifically allowed to do so in this or another Medtronic document. As with all orthopaedic and neurosurgical implants, none of the CD HORIZON™ Spinal System components should ever be reused under any circumstances.

INDICATIONS

The CD HORIZON™ Spinal System with or without SEXTANT™ instrumentation is intended for posterior, non-cervical fixation as an adjunct to fusion for the following indications: degenerative disc disease (DDD - defined as back pain of discogenic origin with degeneration of the disc confirmed by history and radiographic studies), spondylolisthesis, trauma (i.e. fracture or dislocation), spinal stenosis, curvatures (i.e. scoliosis, kyphosis, or lordosis), tumor, pseudarthrosis, and/or failed previous fusion.

Except for hooks, when used as an anterolateral thoracic/lumbar system, the CD HORIZON™ Spinal System may also be used for the same indications as an adjunct to fusion.

With the exception of DDD, the CD HORIZON™ LEGACY™ 3.5mm rods and the CD HORIZON™ Spinal System PEEK rods and associated components may be used for the aforementioned indications in skeletally mature patients as an adjunct to fusion. The 3.5mm rods may be used for the specific pediatric indications noted below.

When used for posterior non-cervical pedicle screw fixation in pediatric patients, the CD HORIZON™ Spinal System implants are indicated as an adjunct to fusion to treat progressive spinal deformities (i.e. scoliosis, kyphosis, or lordosis) including idiopathic scoliosis, neuromuscular scoliosis, and congenital scoliosis. Additionally, the CD HORIZON™ Spinal System is intended to treat pediatric patients diagnosed with the following conditions: spondylolisthesis/spondylolysis, fracture caused by tumor and/or trauma, pseudarthrosis, and/or failed previous fusion. These devices are to be used with autograft and/or allograft. Pediatric pedicle screw fixation is limited to a posterior approach.

The CD HORIZON™ SPIRE™ Plate is a posterior, single-level, non-pedicle supplemental fixation device intended for use in the non-cervical spine (T1-S1) as an adjunct to fusion in skeletally mature patients. It is intended for plate fixation/attachment to spinous processes for the purpose of achieving supplemental fixation in the following conditions: DDD (as previously defined), spondylolisthesis, trauma, and/or tumor.

In order to achieve additional levels of fixation, the CD HORIZON™ Spinal System rods may be connected to the VERTEX™ Reconstruction System with the VERTEX™ rod connector. Refer to the VERTEX™ Reconstruction System Package Insert for a list of the VERTEX™ indications of use.

CONTRAINDICATIONS

Contraindications include:

- Active infectious process or significant risk of infection (immunocompromise).
- Signs of local inflammation.
- Fever or leukocytosis.
- Morbid obesity.
- Pregnancy.
- Mental illness.
- Grossly distorted anatomy caused by congenital abnormalities.
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery such as the presence of congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count.
- Suspected or documented metal allergy or intolerance.
- Any case not needing a bone graft and fusion.
- Any case where the implant components selected for use would be too large or too small to achieve a successful result.
- Any patient having inadequate tissue coverage over the operative site or inadequate bone stock or quality.
- Any patient in which implant utilization would interfere with anatomical structures or expected physiological performance.
- The CD HORIZON™ SPIRE™ Plate and the CD HORIZON™ PEEK Rods are specifically contraindicated for use in pediatric patients.
- Any patient unwilling to follow postoperative instructions.
- Any case not described in the indications.

Nota bene: although not absolute contraindications, conditions to be considered as potential factors for not using this device include:

- Severe bone resorption.
- Osteomalacia.
- Severe osteoporosis.

POTENTIAL ADVERSE EVENTS

All of the possible adverse events associated with spinal fusion surgery without instrumentation are possible. With instrumentation, a listing of potential adverse events includes:

- Early or late loosening of any or all of the components.
- Disassembly, bending, or breakage of any or all of the components.
- Foreign body (allergic) reaction to implants, debris, corrosion products (from crevice, fretting, or general corrosion) including metallosis, staining, tumor formation, or autoimmune disease.
- Pressure on the skin from component parts in patients with inadequate tissue coverage over the implant possibly causing skin penetration, irritation, fibrosis, necrosis, or pain.
- Bursitis.
- Tissue or nerve damage caused by improper positioning and placement of implants or instruments.
- Post-operative change in spinal curvature, loss of correction, height, or reduction.
- Infection.
- Dural tears, pseudomeningocele, fistula, persistent CSF leakage, meningitis.
- Loss of neurological function (e.g. sensory or motor) including paralysis (complete or incomplete), dysesthesias, hyperesthesia, anesthesia, paresthesia, appearance of radiculopathy, or the development or continuation of pain, numbness, neuroma, spasms, sensory loss, tingling sensation, or visual deficits.
- Cauda equina syndrome, neuropathy, neurological deficits (transient or permanent), paraplegia, paraparesis, reflex deficits, irritation, arachnoiditis, muscle loss.
- Urinary retention or loss of bladder control or other types of urological system compromise.
- Scar formation possibly causing neurological compromise or compression around nerves or pain.
- Fracture, microfracture, resorption, damage, or penetration of any spinal bone (including the sacrum, pedicles, or vertebral body) or bone graft or bone graft harvest site at, above, or below the level of surgery.
- Retropulsed graft.
- Herniated nucleus pulposus, disc disruption, or degeneration at, above, or below the level of surgery.
- Non-union (or pseudarthrosis), delayed union, or mal-union.
- Cessation of any potential growth of the operated portion of the spine.
- Loss of or increase in spinal mobility or function.
- Inability to perform the activities of daily living.
- Bone loss or decrease in bone density, possibly caused by stresses shielding.
- Graft donor site complications including pain, fracture, or wound healing problems.
- Ileus, gastritis, bowel obstruction or loss of bowel control, or other types of gastrointestinal system compromise.
- Hemorrhage, hematoma, occlusion, seroma, edema, hypertension, embolism, stroke, excessive bleeding, phlebitis, wound necrosis, wound dehiscence, damage to blood vessels, or other types of cardiovascular system compromise.
- Reproductive system compromise including sterility, loss of consortium, and sexual dysfunction.
- Development of respiratory problems (e.g. pulmonary embolism, atelectasis, bronchitis, pneumonia, etc.)
- Change in mental status.
- Death.

Note: additional surgery may be necessary to correct some of these potential adverse events.

ADDITIONAL POTENTIAL ADVERSE EVENTS FOR PEDIATRIC PATIENTS

- Inability to use pedicle screw fixation due to anatomic limitations (pedicle dimensions and/or distorted anatomy).
- Pedicle screw malpositioning, with or without neurological or vascular injury.
- Proximal or distal junctional kyphosis.
- Pancreatitis.

WARNING

The safety and effectiveness of pedicle screw spinal systems have been established only for spinal conditions with significant mechanical instability or deformity requiring fusion with instrumentation. These conditions are significant mechanical instability or deformity of the thoracic, lumbar, and sacral spine secondary to degenerative spondylolisthesis with objective evidence of neurologic impairment, fracture, dislocation, scoliosis, kyphosis, spinal tumor, and failed previous fusion (pseudarthrosis). The safety and effectiveness of this device for any other conditions are unknown. The implants are not prostheses. In the absence of fusion, the instrumentation and/or one or more of its components can be expected to pull out, bend, or fracture as a result of exposure to every day mechanical stresses.

A device that has been implanted should never be reprocessed or reused under any circumstances. Sterile packaged devices should also never be resterilized. Reprocessing or reuse may compromise the structural integrity of these implants and create a risk of contamination of the implants which could result in patient injury, illness, or death.

ADDITIONAL WARNINGS FOR PEDIATRIC PATIENTS

Warning: the safety and effectiveness of this device has not been established for use as part of a growing rod construct. This device is only intended to be used when definitive fusion is being performed at all instrumented levels.

The use of pedicle screw fixation in the pediatric population may present additional risks when patients are of smaller stature and skeletally immature. Pediatric patients may have smaller spinal structures (pedicle diameter or length) that may preclude the use of pedicle screws or increase the risk of pedicle screw malpositioning and neurological or vascular injury. Patients not skeletally mature that undergo spinal fusion procedures may have reduced longitudinal spinal growth, or may be at risk for rotational spinal deformities (the "crankshaft phenomenon") due to continued differential growth of the anterior spine.

Other adverse events related to pedicle screw fixation, such as screw or rod bending, breakage, or loosening, may also occur in pediatric patients. Pediatric patients may be at increased risk for device-related injury because of their smaller stature.

ADDITIONAL WARNING FOR THE CD HORIZON™ SPIRE™ SPINOUS PROCESS PLATE

Consider the extent of decompression, as well as the amount of intact bone remaining on the spinous processes, when using the CD HORIZON™ SPIRE™ Plate as the sole supplemental fixation for an interbody fusion procedure.

PRECAUTIONS

The implantation of pedicle screw spinal systems should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.

A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise the results. This device system is not intended to be the sole means of spinal support. Use of this product without a bone graft or in cases that develop into a non-union will not be successful. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly, and/or breakage of the device(s) will eventually occur.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and compliance of the patient will greatly affect the results. Patients who smoke have been shown to have an increased incidence of non-unions. These patients should be advised of this fact and warned of this consequence. Obese, malnourished, or alcohol abuse patients are also poor candidates for spine fusion. Patients with poor muscle and bone quality and/or nerve paralysis are also poor candidates for spine fusion.

ADDITIONAL PRECAUTIONS FOR PEDIATRIC PATIENTS

The implantation of pedicle screw spinal systems in pediatric patients should be performed only by experienced spinal surgeons with specific training in the use of this pedicle screw spinal system in pediatric patients because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful utilization of the system in pediatric patients.

The selection of the proper size, shape, and design of the implant for each patient is crucial to the safe use of this device in pediatric patients.

Physician note: although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

For US audiences only

Caution: Federal law (USA) restricts these devices to sale by or on the order of a physician.

IMPLANT SELECTION

The selection of the proper size, shape, and design of the implant for each patient is crucial to the success of the procedure. Metallic surgical implants are subject to repeated stresses in use, and their strength is limited by the need to adapt the design to the size and shape of human bones. Unless great care is taken in patient selection, proper placement of the implant, and postoperative management to minimize stresses on the implant, such stresses may cause metal fatigue and consequent breakage, bending, or loosening of the device before the healing process is complete which may result in further injury or the need to remove the device prematurely.

DEVICE FIXATION

In cases where a percutaneous posterior approach is used, refer to the CD HORIZON™ SEXTANT™ surgical technique. CD HORIZON™ Spinal System instrumentation contains rods and implants of various diameters, which are intended to be used with device specific instruments. For self breaking plugs, always hold the assembly with the Counter Torque device. Tighten and break-off the head of the plug to leave the assembly at optimum fixation security. After the upper part of the self-breaking plug has been sheared off, further re-tightening is not necessary and not recommended. The head part should not remain in the patient. After the upper part of the self breaking plug has been sheared off, re-adjustment is not possible unless the plug is removed and replaced with a new one.

When using DTT Transverse Links, the M6 plug should be tightened between 8 and 9 Nm. (70 to 80 inch-lbs).

CD HORIZON™ PEEK Rods are not to be used with CROSSLINK™ Plates or in pediatric patients.

PREOPERATIVE

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and pre dispositions such as those addressed in the aforementioned contraindications should be avoided.
- Care should be used in the handling and storage of the implant components. The implants should not be scratched or otherwise damaged. Implants and instruments should be protected during storage, especially from corrosive environments.
- An adequate inventory of implants should be available at the time of surgery; normally a quantity in excess of what is expected to be used.
- Since mechanical parts are involved, the surgeon should be familiar with the various components before using the equipment and should personally assemble the devices to verify all parts and necessary instruments are present before surgery. The CD HORIZON™ Spinal System components (described in the DESCRIPTION section) are not to be combined with the components from another manufacturer.
- Unless sterile packaged, all devices should be sterilized before use. Additional sterile components should be available in case of an unexpected need.

INTRAOPERATIVE

- Extreme caution should be used around the spinal cord and nerve roots. Damage to the nerves will cause loss of neurological functions.
- Breakage, slippage, or misuse of instruments or implant components may cause injury to the patient or operative personnel.
- The rods should not be repeatedly or excessively bent. The rods should not be reverse-bent in the same location. Use great care to ensure the implant surfaces are not scratched or notched since such actions may reduce the functional strength of the construct. If the rods are cut to length, they should be cut in such a way as to create a flat, non-sharp surface perpendicular to the midline of the rod. Cut the rods outside the operative field. Whenever possible, use pre-cut rods of the length needed.
- Utilize an imaging system to facilitate surgery.
- To insert a screw properly, a guide wire should first be used, followed by a sharp tap. Caution: be careful the guide-wire, if used, is not inserted too deep, becomes bent, and/or breaks. Ensure the guide-wire does not advance during tapping or screw insertion. Remove the guide-wire and make sure it is intact. Failure to do so may cause the guide wire or part of it to advance through the bone and into a location that may cause damage to underlying structures.
- Caution: do not overlap or use a screw/bolt that is either too long or too large. Overlapping, using an incorrectly sized screw/bolt, or accidentally advancing the guidewire during tap or screw/bolt insertion may cause nerve damage, hemorrhage, or the other possible adverse events listed elsewhere in this package insert. If screws/bolts are being inserted into spinal pedicles, use as large a screw/bolt diameter as will fit into each pedicle.
- Bone graft must be placed in the area to be fused and graft material must extend from the upper to the lower vertebrae being fused.
- To ensure maximum stability, two or more CROSSLINK™ plates or DTT Transverse Links on two bilaterally placed, continuous rods should be used whenever possible.
- Bone cement should not be used because the safety and effectiveness of bone cement has not been determined for spinal uses, and this material will make removal of the components difficult or impossible. The heat generated from the curing process may also cause neurologic damage and bone necrosis.
- Before closing the soft tissues, provisionally tighten (finger tighten) all of the nuts or screws, especially screws or nuts that have a break-off feature. Once this is completed, go back and firmly tighten all of the screws and nuts. Recheck the tightness of all nuts or screws after finishing to make sure none loosened during the tightening of the other nuts or screws. Failure to do so may cause loosening of the other components.

POSTOPERATIVE

The physician's postoperative directions and warnings to the patient, and the corresponding patient compliance, are extremely important.

- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening, and/or breakage of the device(s) are complications which may occur as a result of excessive or early weight-bearing or muscular activity. The risk of bending, loosening, or breakage of a temporary internal fixation device during postoperative rehabilitation may be increased if the patient is active, or if the patient is debilitated or demented. The patient should be warned to avoid falls or sudden jolts in spinal position.
- To allow the maximum chances for a successful surgical result, the patient or devices should not be exposed to mechanical vibrations or shock that may loosen the device construct. The patient should be warned of this possibility and instructed to limit and restrict physical activities, especially lifting and twisting motions and any type of sport participation. The patient should be advised not to smoke tobacco, utilize nicotine products, or consume alcohol or non-steroidals or anti-inflammatory medications such as aspirin during the bone graft healing process.
- The patient should be advised of their inability to bend or rotate at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
- Failure to immobilize a delayed or non-union of bone will result in excessive and repeated stresses on the implant. By the mechanism of fatigue, these stresses can cause the eventual bending, loosening, or breakage of the device(s). It is important that immobilization of the spinal surgical site be maintained until firm bony union is established and confirmed by roentgenographic examination. If a state of non-union persists or if the components loosen, bend, or break, the device(s) should be revised or removed immediately before serious injury occurs. The patient must be adequately warned of these hazards and closely supervised to ensure cooperation until bony union is confirmed.
- As a precaution, before patients with implants receive any subsequent surgery (such as dental procedures), prophylactic antibiotics may be considered, especially for high-risk patients.
- The CD HORIZON™ Spinal System implants are temporary internal fixation devices. Internal fixation devices are designed to stabilize the operative site during the normal healing process. After the spine is fused, these devices serve no functional purpose and may be removed. While the final decision on implant removal is, of course, up to the surgeon and patient, in most patients removal is indicated because the implants are not intended to transfer or support forces developed during normal activities. If the device is not removed following completion of its intended use, one or more of the following complications may occur: (1) corrosion, with localized tissue reaction or pain; (2) migration of implant position possibly resulting in injury; (3) risk of additional injury from postoperative trauma; (4) bending, loosening, and breakage which could make removal impractical or difficult; (5) pain, discomfort, or abnormal sensations due to the presence of the device; (6) possible increased risk of infection; (7) bone loss due to stress shielding; and (8) potential unknown or unexpected long term effects such as carcinogenesis. Implant removal should be followed by adequate postoperative management to avoid fracture, re-fracture, or other complications.
- Any retrieved devices should be treated in such a manner that reuse in another surgical procedure is not possible. As with all orthopedic implants, the CD HORIZON™ Spinal System components should never be reused under any circumstances.

PACKAGING

Devices may be supplied in a sterile or non-sterile form. Packages for each of the components should be intact upon receipt. Once the seal on the sterile package has been broken, the product should not be re-sterilized. If a loaner or consignment system is used, all sets should be carefully checked for completeness and all components should be carefully checked to ensure there is no damage prior to use. Damaged packages or products should not be used, and should be returned to Medtronic.

STERILIZATION

Unless marked sterile and clearly labeled as such in an unopened sterile package provided by the company, all implants used in surgery must be sterilized by the hospital prior to use. Remove all packaging materials prior to sterilization. Only sterile products should be placed in the operative field. Unless specified elsewhere, these products are recommended to be steam-sterilized by the hospital using one of the sets of process parameters in Table 1.

Table 1: Sterilization cycle parameters for the US and its territories

Method	Cycle	Temperature	Exposure time	Minimum dry time ¹
Steam	Dynamic-air-removal	270°F (132°C)	4 Minutes	30 Minutes
Steam	Dynamic-air-removal	275°F (135°C)	3 Minutes	30 Minutes

Note: because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g. temperatures, times) used for their equipment. It is the end user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) that have been cleared by the Food and Drug Administration (FDA) for the selected sterilization cycle specifications (time and temperature). The sterilization cycles listed in Table 2 are not considered by the FDA to be standard sterilization cycles.

*For medical facilities located outside the US and its territories, some non-US health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical instruments that could come into contact with the central nervous system.

Table 2: Sterilization cycle parameters for medical facilities outside the US and its territories

Method	Cycle	Temperature	Exposure time	Minimum dry time ¹
Steam	Dynamic-air-removal	273°F (134°C)	4 Minutes	30 Minutes
Steam	Dynamic-air-removal	273°F (134°C)	20 Minutes*	30 Minutes

¹The minimum dry times were validated using sterilizers having vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer's recommendations.

Only sterile products should be placed in the operative field. The general instruments used with this device are provided non-sterile. Refer to the instrument package insert for sterilization parameters and requirements. No implant should be re-used once it comes into contact with human tissue or body fluid. Always immediately clean and re-sterilize instruments that have been used in surgery. This process must be performed before handling or (if applicable) returning to Medtronic.

MRI INFORMATION

The CD HORIZON™ Spinal System has not been evaluated for safety, heating, migration, or compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the CD HORIZON™ Spinal System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

PRODUCT COMPLAINTS

For product problems, contact Medtronic.

FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is needed or required, contact Medtronic.

©2018 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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IMPORTANT INFORMATION ON THE MEDTRONIC NAVIGATED REUSABLE INSTRUMENTS FOR USE WITH STEALTHSTATION™ AND IPC™ POWEREASE™ SYSTEMS

DESCRIPTION

Medtronic Navigated Reusable Instruments are spine preparation instruments made of high grade stainless steel. These instruments are specifically designed for use in procedures where the use of stereotactic surgery may be appropriate. Placing Medtronic single-use sterile spheres on each of the NavLock™ Tracker passive stems allows a Medtronic computer-assisted surgery system such as the StealthStation™ Image Guidance System to track the instruments in the surgical field.

Medtronic Navigated Reusable Instruments are compatible with various Medtronic spinal implant systems. These instruments are also compatible with Medtronic's IPC™ POWEREASE™ System when connected to the POWEREASE™ Driver.

INTENDED USE

Medtronic Navigated Reusable Instruments are intended to be used during the preparation and placement of Medtronic screws during spinal surgery to assist the surgeon in precisely locating anatomical structures in either open, or minimally invasive, procedures. Medtronic Navigated Reusable Instruments are specifically designed for use with the StealthStation™ System, which is indicated for any medical condition in which the use of stereotactic surgery may be appropriate, and where reference to a rigid anatomical structure, such as a skull, a long bone, or vertebra can be identified relative to a CT or MR-based model, fluoroscopy images, or digitized landmarks of the anatomy. Medtronic Navigated Reusable Instruments are also compatible with the IPC™ POWEREASE™ System.

Do not implant the instruments.

If there is any doubt or uncertainty concerning the proper use of these instruments, contact Medtronic. Any available surgical techniques will be provided at no charge.

Medtronic does not and cannot warrant the use of this instrument nor any of the component parts upon which repairs have been made or attempted, except as performed by Medtronic or an authorized Medtronic repair representative. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

WARNINGS

- Breakage, slippage, misuse, or mishandling of instruments, such as on sharp edges, may cause injury to the patient or operative personnel.
- Improper maintenance, handling, or poor cleaning procedures can render the instrument unsuitable for its intended purpose or even dangerous to the patient or surgical staff.
- It is important the surgeon exercise extreme caution when working in close proximity to vital organs, nerves, or vessels and the forces applied while correcting the position of the instrumentation is not excessive, such that it might cause injury to the patient.
- During navigation, it is important to frequently confirm navigational accuracy by touching the tip of the instrument on known anatomical points, including accuracy checkpoints, and comparing the position of the instrument tip in the image with its physical location.

PRECAUTIONS

- Excessive force applied by instruments to implants can dislodge devices, particularly hooks.
- Never expose instruments to temperatures in excess of 135°C that may modify the physical characteristics.
- Extreme care should be taken to ensure that this instrument remains in good working order. During the procedure, successful utilization of this instrument is extremely important. Instruments should not be bent or damaged in any way. Misuse of instruments resulting in corrosion, "freezing-up", scratching, loosening, bending, or fracture of any or all sections of an instrument may inhibit or prevent proper function.
- These instruments should be carefully placed on trays, cleaned after each use, and stored in a dry environment.
- Do not use this instrument for any action for which it was not intended.
- Regularly review the operational state of all instruments and, if necessary, make use of repair and replacement services.
- To avoid injury or navigation inaccuracy, the instrument should be carefully examined for functionality or damage prior to use. A damaged instrument should not be used. Additional back-up instruments should be available.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful utilization of the system by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results.
- Proper patient selection and operative care are critical to the success of the surgery and avoidance of injury during surgery. Read and follow all other product information supplied by the manufacturer of the implants or the instruments.
- Special precautions are needed during pediatric use. Care should be taken when using instruments in pediatric patients since these patients can be more susceptible to the stresses involved in their use.
- Some surgeries require the use of instruments which incorporate a measuring function. Ensure these are not worn and any surface engravings are clearly visible.

POTENTIAL ADVERSE EFFECTS

- Nerve damage, paralysis, pain, or damage to soft tissue, visceral organs, or joints.
- Infection if instruments are not properly cleaned and sterilized.
- Pain, discomfort, or abnormal sensations resulting from the presence of the instrument.
- Nerve damage due to surgical trauma.
- Dural leak in cases of excessive load application.
- Impingement or damage of close vessels, nerves, and organs by slippage or misplacement of the instrument.
- Damage due to spontaneous release of clamping devices or spring mechanisms of certain instruments.
- Cutting of skin or gloves of operating staff.
- Bony fracture in cases of deformed spine or weak bone.
- Tissue damage to the patient, physical injury to operating staff, and/or increased operating time that may result from the disassembly of multi-component instruments occurring during surgery.
- The methods of use of instruments are determined by the user's experience and training in surgical procedures. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results.

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This device should be used only by physicians familiar with the device, its intended use, any additional instrumentation, and any available surgical techniques.

PACKAGING

Packages for components should be intact upon receipt. All sets should be checked for completeness and all components should be checked for signs of damage prior to use. Damaged packages or products should not be used and should be returned to Medtronic.

All instruments used in surgery must be cleaned and sterilized by the hospital prior to use.

Remove all packaging material prior to sterilization. Only sterile instruments should be used in surgery. Always immediately clean and re-sterilize all instruments used in surgery. Instruments should be thoroughly cleaned prior to re-sterilization. This process must be performed before handling, or before returning product to Medtronic.

EXAMINATION

Instruments must always be examined by the user prior to surgery.

Examination should be thorough and must include a visual and functional inspection of the working surfaces, pivots, racks, spring or torsional operation and the presence of any cracks, bending, deformation, or distortion, and that all components are complete.

Never use instruments with obvious signs of damage or that are incomplete or otherwise nonfunctional.

Visual Inspection

Make certain of the following:

- Laser etchings, engravings, and other markings are legible.
- No cracks are present in instrument handles or any part of the instrument.
- Discoloration, corrosion, stains, or rust do not exist.
- There is no handle/shaft separation, and the handle-to-shaft connection is secure.
- No cuts or gouges in silicone are present.
- There is no damage (cuts, tears, etc.) to insulation.
- There is no damage to the working ends or tips. The working ends should be free of cracks, sharp edged gouges, and other damage. When applicable, the working ends should be sharp.
- The instrument tip and/or shaft is not bent.
- There is no damage to threads.
- All parts are present and free of damage and deterioration. Examples of parts that may be missing, loose, or damaged include set screws, springs, curved springs, pins, and prongs.
- Mating ends are free of damage (nicks, gouges, bends, etc.) that would interfere with the mating function.
- Cannulated instruments with a guide wire or other insertion tool are visually checked.

Functional Inspection

Make certain of the following:

- The parts intended to move do so freely without sticking, binding, or grinding.
- Springs return the handle of the instrument to its original position.
- Retention tabs hold appropriate mating parts and are not damaged.
- The instrument will function as intended with the appropriate mating parts.
- Ball detents will hold mating parts and are free from damage.
- Sharp edges are sharp to the touch and are not dull, have no nicks, or any other damage.
- Tips meet when appropriate.
- Ratcheting mechanisms are functional. This includes handles, latches, and other mechanisms. All teeth should be present and functional.
- Driver tips are not worn beyond functional use. If necessary, mate the instrument with the appropriate part.

DIRECTIONS FOR USE

The instruments are precision devices, which may incorporate a measuring function and have uses as described on the label. Those instruments which incorporate a depth measuring function have an accuracy level equivalent to $\pm 0.1\text{mm}$.

Unless labeled for single use, the instruments may be re-used.

Use with the IPC™ POWEREASE™ System

See the IPC™ POWEREASE™ System package insert (M726750B246) for instructions on how to set up and operate the IPC™ POWEREASE™ System.

Using the POWEREASE™ Adapter is optional. See the POWEREASE™ Adapter package insert (9734925) for instructions on using the adapter with the POWEREASE™ Driver.

1. Place single-use sterile spheres on each of the four NavLock™ Tracker stems. Push each sphere onto the stem until it "clicks" into place. Make sure the sphere is firmly seated on the stem.
2. Place the instrument into the NavLock™ Tracker.
 - A. Insert the proximal end of the instrument shaft until it is fully seated in the NavLock™ Tracker.
 - B. Verify the tool is secured in the NavLock™ Tracker by making sure it cannot be pulled out of the device.
3. Attach the NavLock™ Tracker assembly to the POWEREASE™ Driver (Figure 1).
 - A. On the POWEREASE™ Driver, pull back and hold the quick disconnect to unlock the collet.

Note: The quick disconnect must be held in the unlocked position ① when inserting or removing tools.

- B. Insert the proximal end of the NavLock™ Tracker assembly into the collet.
- C. Align the flat sides of the instrument shaft with the marks on the collet ⑧ and insert into the collet until fully seated (Figure 2).
- D. Release the quick disconnect to place the assembly in the locked position ②.
- E. Verify the NavLock™ Tracker assembly is secured in the POWEREASE™ Driver by ensuring it cannot be pulled out of the device.

Table 1: POWEREASE™ Driver

①	Quick Disconnect Unlocked	⑤	Finger Sensor Trigger
②	Quick Disconnect Locked	⑥	NIM-ECLIPSE Cable Connector
③	Mechanized Working Collar Locking Pins (2)	⑦	Mode Select Switch
④	Variable Speed Trigger	⑧	Flat Side Indicators

- F. Refer to the Synergy™ Spine & Trauma Pocket Guide for additional navigation instructions.
4. Verify the assembly.
 - A. For Synergy™ Spine Version 1.7, select the current NavLock™ Tracker color and corresponding instrument from the instrument list (tool card) in the procedure software. For Synergy™ Spine & Trauma software version 2.0 and later, select the appropriate NavLock™ tool card.
 - B. Ensure the instrument has green status (i.e., it is being tracked by the StealthStation™ system) on the screen.
 - C. Face the NavLock™ Tracker array toward the camera.
 - D. Place the instrument tip straight into the reference frame divot.
 - E. Press the footswitch or wait for the system to verify.
5. In the procedure software, select the appropriate instrument that corresponds to the instrument being used in the procedure.
 - A. While in the Navigate task, click the Select Tip button.
 - B. Select the instrument type and size from the list in the Select Tip flyout.
6. Refer to the Synergy™ Spine & Trauma Pocket Guide for additional navigation instructions.

When using the POWEREASE™ Driver with navigation, avoid any bending of the instrument by applying direct force only along the instrument assembly axis (Figure 3).

Instructions for Manual Use

1. Place single-use sterile spheres on each of the four NavLock™ Tracker stems. Push each sphere onto the stem until it "clicks" into place. Ensure the sphere is firmly seated on the stem.
2. Place the instrument into the NavLock™ Tracker.
 - A. Insert the proximal end of the instrument shaft until it is fully seated in the NavLock™ Tracker.
 - B. Verify the instrument is secured in the NavLock™ Tracker by ensuring it cannot be pulled out of the device.
3. Secure the handle to the NavLock™ Tracker assembly. Note: Use only the Medtronic Ratcheting Egg Handle (9734410), Medtronic Straight Handle (9733734), Quick Connect Ratcheting Handle (G900000), or Medtronic Ratcheting Handle (G170059).
 - A. Retract the collar on the handle and snap the proximal end of the instrument shaft into the handle.
 - B. Verify the handle is secured to the instrument by ensuring it cannot be pulled out of the device.
4. Continue with steps 4 through 6 in the section entitled "Use with the IPC™ POWEREASE™ System."

REPROCESSING - GENERAL CONSIDERATIONS

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- Minimizes the organic soil transfer from one patient to another.
- Prevents accumulation of residual soil throughout the product's use life.
- Allows for successful follow up sterilization. Adequate reprocessing is contingent upon the thoroughness of cleaning.

Cleaning is the initial step and sterilization occurs later in reprocessing and is intended to kill microorganisms to reduce the likelihood of transmission and possibilities of infection. To ensure acceptable reprocessing, there should be no delay between the steps in this document.

Bloodborne Pathogens

Universal precautions for handling this instrument after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030 regarding occupational exposure to bloodborne pathogens.

Thorough Cleaning of Instruments

It is critical that all reusable instruments are thoroughly cleaned by the hospital prior to use and after each use following the specific cleaning steps listed in this document. Ineffective or incomplete cleaning can hinder subsequent sterilization activities.

Automated Cleaning

A fully automated washer/disinfector cycle is not recommended as the sole cleaning method. Manual cleaning is required.

Manual Cleaning

Disinfection agents such as sanitizing and chemical solutions act to reduce microbes on surfaces, but may not act as cleaners. Other types of soaps or detergents may not clean to an appropriate level of cleanliness. Enzymatic cleaning agents are recommended. This type of cleaner has been shown to effectively remove organic soils, such as blood, from instruments. The cleaning instructions have been validated by Medtronic using an enzymatic cleaner.

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The following cleaning agents, solutions, or tools should not be used:

- Saline solution.
- Alkaline cleaning agents.
- Solutions containing chlorine (e.g., bleach) or aldehydes (e.g., glutaraldehyde).
- Formalin, mercury, chlorides, bromides, iodides, or Ringer's solution.
- Metal brushes or scouring pads.

The use of neutral pH enzymatic cleaners and soft bristled brushes and soft pipe cleaners are recommended.

Proper Handling After Use

Do not allow the instruments to dry after use and prior to cleaning. Cleaning and subsequent sterilization may be hindered when blood or bloody solutions are allowed to dry on instruments.

Cleaning Instructions - Point of Use

1. Remove all visible soil from instruments using non-shedding wipes.
2. Place the instruments in a tray of water or cover with damp towels. Instruments should be cleaned within 30 minutes of use to minimize the potential for drying.

Cleaning Instructions - Dedicated Cleaning Area

1. Immediately transport the tray containing covered instruments to a work area dedicated to further reprocessing.
2. Rinse instruments under running tap water for 3 minutes.
3. Scrub instruments with an appropriate soft bristle brush to aid in the removal of visible soil. Scrub inside any lumens or cavities.
4. Repeat steps 2 and 3 until all visible soil is removed.
5. Using tap water, prepare an enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures.
6. Place instruments in the enzymatic cleaner, completely submerged, and soak for 45-60 minutes.
7. Remove instruments from the enzymatic cleaner and rinse under running tap water for 3 minutes. Flush any lumens or cavities in the water stream.
8. Using tap water, prepare a second enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures in an appropriately-sized sonicator.
9. Place instruments in the enzymatic cleaner, completely submerged, and sonicate for 45-60 minutes.
10. Remove instruments from the sonicator and rinse using running tap water for 3 minutes. Flush any lumens or cavities in the water stream.
11. Repeat rinsing step 10, using de-ionized water for an additional 3 minutes.
12. Dry the instruments using clean, absorbent, non-shedding wipes.
13. Carefully inspect instruments, including any lumens and cavities, to ensure all contamination has been removed. If any soil is still present, repeat the cleaning process or contact your Medtronic distributor immediately to arrange for disposal or replacement. Do not proceed with reprocessing of a soiled instrument.

Disinfection

The Working Ends must be sterilized prior to initial use, or in adherence to these reprocessing instructions before re-use. Disinfection using solutions or chemicals is unnecessary and not recommended.

Sterilization Instructions

Only sterile products should be placed in the operative field. Unless marked sterile and clearly labeled as such in an unopened sterile package provided by the company, all instruments used in surgery must be sterilized by the hospital prior to use.

Ethylene oxide (EO), gas plasma, gamma irradiation, chemical vapor, or dry heat sterilization methods are NOT recommended for sterilization of these instruments. Steam/moist heat is the recommended method of sterilization.

Devices must be sterilized prior to initial use, or as part of these reprocessing instructions, before re-use.

The following sterilization instructions have been validated to a sterility assurance level of 10^{-6} .

1. Wrap or pouch the instruments.
2. Inspect the packaging to ensure no rips, punctures, or seal failures are present in or on the packaging prior to loading into the sterilizer.
3. Load the instruments into the sterilizer by following the sterilizer manufacturer's recommended loading procedures and load configurations.
4. Follow the sterilizer manufacturer's recommended procedures to program the sterilizer with any one of the sets of sterilization cycle parameters found in Table 2 and Table 3.

Table 2: Sterilization cycle parameters for the US and its territories

Method	Cycle	Temperature	Exposure time	Minimum dry time ¹
Steam	Gravity Displacement	250°F (121°C)	30 Minutes	30 Minutes
Steam	Gravity Displacement	270°F (132°C)	15 Minutes	30 Minutes
Steam	Gravity Displacement	275°F (135°C)	10 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	270°F (132°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	275°F (135°C)	3 Minutes	30 Minutes

For medical facilities located outside the US and its territories: Some non-US health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical instruments that could come into contact with the central nervous system.

Table 3: Sterilization cycle parameters for medical facilities outside the US and its territories

Method	Cycle	Temperature	Exposure time	Minimum dry time ¹
Steam	Gravity Displacement	273°F (134°C)	20 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	20 Minutes	30 Minutes

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NOTE: Chamber size and chamber load differences may exist between industrial and health care facility sterilizer models. The sterilization parameters listed in Tables 2 and 3 can be achieved in both health care facility and larger, industrial sterilizer models. Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g., temperatures, times) used for their equipment.

The sterilization cycles listed in Table 3 are not considered by the Food and Drug Administration (FDA) to be standard sterilization cycles. It is the user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, chemical indicators, biological indicators, and sterilization cassettes) cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

FURTHER INFORMATION

In case of complaint, or for supplementary information, contact Medtronic.

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Medtronic Reusable Instruments for use with the IPC® POWEREASE® System

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2015-05-22

AUSTRALIAN SPONSOR:

Medtronic Australia Pty Ltd
97 Waterloo Rd
North Ryde, NSW 2113
Australia

Medtronic Sofamor Danek USA, Inc.
1000 Pyramid Place
Memphis, TN 38132
Telephone: 800 933 2635 (USA)
901 396 3133 (Outside USA)
Fax: 901 396 8356

Medtronic B.V.
Earl Bakkenstraat 10
6422 PZ Heerlen
The Netherlands
Tel: + 31 45 566 80 00

ENGLISH

IMPORTANT INFORMATION FOR MEDTRONIC REUSABLE INSTRUMENTS FOR USE WITH THE IPC® POWEREASE® SYSTEM

DESCRIPTION

The Medtronic Reusable Instruments compatible with Medtronic's IPC® POWEREASE® System are spine preparation instruments, manufactured from materials commonly used in orthopedic procedures which meet available national or international standards specifications. The instruments may be connected to the POWEREASE® Driver or used manually. These instruments are also compatible with various Medtronic spinal implant systems.

INDICATIONS FOR USE

IPC® System is indicated for the incision/cutting, removal, drilling and sawing of soft and hard tissue and bone, and biomaterials in Neurosurgical (Cranial, Craniofacial), Orthopedic, Arthroscopic, Spinal, Stereotaxy, and General surgical procedures.

The IPC® POWEREASE® System is indicated for drilling, tapping and driving screws and working end attachments during spinal surgery, including open and minimally invasive procedures. It is also used in placement or cutting of screws, posts and rods.

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Medtronic does not and cannot warrant the use of this instrument, nor any of the instrument component parts upon which repairs have been made or attempted, except as performed by Medtronic or an authorized Medtronic repair representative. Implied warranties of merchantability and fitness for a particular purpose or use are specifically excluded.

WARNINGS

- Breakage, slippage, misuse, or mishandling of instruments, such as on sharp edges, may cause injury to the patient or operative personnel.
- Improper maintenance, handling, or poor cleaning procedures can render the instrument unsuitable for its intended purpose or even dangerous to the patient or surgical staff.
- It is important the surgeon exercise extreme caution when working in close proximity to vital organs, nerves, or vessels and the forces applied while correcting the position of the instrumentation is not excessive, such that it might cause injury to the patient.

PRECAUTIONS

- Excessive force applied by instruments to implants can dislodge devices.
- Never expose instruments to temperatures in excess of 135°C that may considerably modify the physical characteristics.
- Extreme care should be taken to ensure these instruments remain in good working order. During the procedure, successful utilization of these instruments is extremely important. Instruments should not be bent or damaged in any way. Misuse of instruments resulting in corrosion, "freezing-up", scratching, loosening, bending, or fracture of any or all sections of an instrument may inhibit or prevent proper function.
- These instruments should be carefully placed on trays, cleaned after each use, and stored in a dry environment.
- Do not use this instrument for any action for which it was not intended.
- Regularly review the operational state of all instruments and, if necessary, make use of repair and replacement services.
- To avoid injury, the instrument should be carefully examined for functionality or damage prior to use. A damaged instrument should not be used. Additional back-up instruments should be available.
- Preoperative and operating procedures, including knowledge of surgical techniques, are important considerations in the successful utilization by the surgeon. Further, the proper selection and the compliance of the patient will greatly affect the results.
- Proper patient selection and operative care are critical to the success of the surgery and avoidance of injury during surgery. Read and follow all other product information supplied by the manufacturer of the implants or the instruments.
- Special precautions are needed during pediatric use. Care should be taken when using instruments in pediatric patients since these patients can be more susceptible to the stresses involved in their use.
- Some surgeries require the use of instruments which incorporate a measuring function. Ensure that these are not worn and any surface engravings are clearly visible.

POSSIBLE ADVERSE EFFECTS

- Nerve damage, paralysis, pain, or damage to soft tissue, visceral organs, or joints.
- Infection if instruments are not properly cleaned and sterilized.
- Pain, discomfort, or abnormal sensations resulting from the presence of the instrument.
- Nerve damage due to surgical trauma.
- Dural leak in cases of excessive load application.
- Impingement of close vessels, nerves, and organs by slippage or misplacement of the instrument.
- Damage due to spontaneous release of clamping devices or spring mechanisms of certain instruments.
- Cutting of skin or gloves of operating staff.
- Bony fracture in cases of deformed spine or weak bone.
- Tissue damage to the patient, physical injury to operating staff, and/or increased operating time that may result from the disassembly of multi-component instruments occurring during surgery.
- The methods of use of instruments are determined by the user's experience and training in surgical procedures. A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where other patient conditions may compromise the results.

PHYSICIAN NOTE: Although the physician is the learned intermediary between the company and the patient, the important medical information given in this document should be conveyed to the patient.

USA For US audiences only

CAUTION: FEDERAL LAW (USA) RESTRICTS THESE DEVICES TO SALE BY OR ON THE ORDER OF A PHYSICIAN.

These instruments should be used only by physicians familiar with the device, its intended use, any additional instrumentation, and any available surgical techniques.

PACKAGING

Packages for instrument components should be intact upon receipt. All sets should be checked for completeness and all instrument components should be checked for signs of damage prior to use. Damaged packages or products should not be used and should be returned to Medtronic.

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Remove all packaging material prior to cleaning and sterilization. Only sterile instruments should be used in surgery. Always immediately clean and re-sterilize all instruments used in surgery. Instruments should be thoroughly cleaned prior to re-sterilization. This process must be performed before handling, or before returning product to Medtronic.

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Examination should be thorough and must include a visual and functional inspection of the working surfaces, spring or torsional operation, cleanliness of location holes or cannulations, and the presence of any cracks, bending, deformation, or distortion, and that all components are complete.

Never use instruments with obvious signs of excessive wear, damage, or that are incomplete or otherwise nonfunctional.

Visual Inspection

Make certain of the following:

- Laser etchings, engravings, and other markings are legible.
- No cracks are present in instrument handles or any part of the instrument.
- Discoloration, corrosion, stains, or rust do not exist. If present, attempt to remove in accordance with the cleaning instructions provided in the Reprocessing section of this document.
- The handle-to-shaft connection is secure.
- No cuts or gouges in silicone are present.
- There is no damage (cuts, tears, etc.) to insulation.
- There is no damage to the working end of the instrument. The instrument should be free of cracks, sharp edged gouges, and other damage. When applicable, the working ends should be sharp.
- There is no damage to threads.
- All parts are present and free of damage and deterioration.
- Mating ends are free of damage (nicks, gouges, bends, etc.) that would interfere with the mating function.
- Cannulated instruments with a guide wire or other insertion tool are visually checked.

Functional Inspection

Make certain of the following:

- The parts intended to move will do so freely without sticking, binding, or grinding.
- Retention mechanism holds appropriate mating parts and is not damaged.
- The instrument will function as intended with the appropriate mating parts.
- Ball detents will hold mating parts and are free from damage.
- Sharp edges are sharp to the touch and are not dull, have no nicks, or any other damage.
- Driver tips are not worn beyond functional use. If necessary, mate the instrument with the appropriate part.

DIRECTIONS FOR USE

The instruments are precision devices which may incorporate a measuring function and have uses as described on the label.

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Use with the IPC® POWEREASE® System.

See the IPC® POWEREASE® System package insert (M726750B246) for instructions on setting up and operating the system.

1. Attach the reusable instrument to the POWEREASE® Driver (Figure 1).
 - a. On the POWEREASE® Driver, pull back and hold the quick disconnect to unlock the collet.
 - Note:** The quick disconnect must be held in the unlocked position ① when inserting or removing tools.
 - b. Insert the proximal end of the instrument into the collet.
 - c. Align the flat sides of the instrument shaft with the marks on the collet ② and insert into the collet until fully seated (Figure 2).
 - d. Release the quick disconnect to place the assembly in the locked position ③. Verify the instrument is secured in the POWEREASE® Driver by ensuring it cannot be pulled out of the device.

When using the POWEREASE® Driver with reusable instruments, avoid any bending of the instrument by applying direct force only along the instrument assembly axis (Figure 3).

Instructions for Manual Use

1. Secure the manual handle to the connection end of the instrument: Use only the Medtronic Quick Connect Ratcheting Handle (G900000) or Medtronic Quick Connect T-Handle (G900100).
2. Retract the collar on the handle and snap the proximal end of the instrument shaft into the handle.
3. Verify the handle is secured to the instrument by ensuring it cannot be pulled out of the device.

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Bloodborne Pathogens

Universal precautions for handling this instrument after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030, Occupational Exposure to Bloodborne Pathogens.

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The following cleaning agents, solutions, or tools should **not** be used:

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Cleaning and Rinsing Water

If available, softened tap water should be used. De-ionized water should be used for the final rinse step to prevent mineral deposits on instrument surfaces.

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These types of lubricants should never be used as they may not be removed by these cleaning instructions. The lubricants may coat microorganisms, prevent direct contact of steam with instrument surfaces, and hinder sterilization.

MANUAL CLEANING INSTRUCTIONS

Proper Handling After Use

Do **not** allow the instruments to dry after use and prior to cleaning. Cleaning and subsequent sterilization may be hindered when blood or bloody solutions are allowed to dry on instruments.

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1. Immediately transport the tray containing covered instruments to a work area dedicated to further reprocessing.
2. Rinse instruments under running tap water for 3 minutes.
3. Scrub instruments with an appropriate soft bristle brush to aid in the removal of visible soil.
4. Repeat steps 2 and 3 until all visible soil is removed.
5. Using tap water, prepare an enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures.
6. Place instruments in the enzymatic cleaner, completely submerged, and soak for 45-60 minutes.
7. Remove instruments from the enzymatic cleaner and rinse under running tap water for 3 minutes. Ensure lumens and cavities are flushed in the water stream.
8. Using tap water, prepare a second enzymatic cleaning solution according to the manufacturer's instructions, dilution recommendations, and temperatures in an appropriately-sized sonicator.
9. Place instruments in the enzymatic cleaner, completely submerged, and sonicate for 45-60 minutes.
10. Remove instruments from the sonicator and rinse using running tap water for 3 minutes. Ensure lumens and cavities are flushed in the water stream.
11. Repeat step 10, using de-ionized water for an additional 3 minutes.
12. Dry the instruments using clean, absorbent, non-shedding wipes.
13. Carefully inspect instruments, including lumens and cavities, to ensure all contamination has been removed. If any soil is still present, repeat the cleaning process or contact your Medtronic distributor immediately to arrange for disposal or replacement. Do not proceed with reprocessing of a soiled instrument.

Inspection

1. Inspect all instruments prior to use. Carefully check for signs of corrosion, discoloring, pitting, or any other signs of wear.
2. Ensure the instruments do not appear bent or damaged. Also inspect the handles for discoloration, cuts, or other signs of damage.

SANITIZATION

Medtronic reusable instruments are considered critical devices and must be cleaned and sterilized prior to initial use, or in adherence to these reprocessing instructions before re-use. Sanitization with disinfectant solutions or chemicals is unnecessary and not recommended.

STERILIZATION INSTRUCTIONS

Only sterile products should be placed in the operative field. Unless marked sterile and clearly labeled as such in an unopened sterile package provided by the company, all instruments used in surgery must be cleaned and sterilized by the hospital prior to use.

Ethylene oxide (EO), gas plasma, gamma irradiation, chemical vapor, or dry heat sterilization methods are **not** recommended for sterilization. Steam/moist heat is the recommended method of sterilization.

Instruments must be cleaned and sterilized prior to initial use, or as part of these reprocessing instructions, before re-use.

It is the user's responsibility to use only sterilizers and accessories (such as sterilization wraps, pouches, chemical indicators, biological indicators, and sterilization trays) cleared by the Food and Drug Administration (FDA) for the selected sterilization cycle parameters (time and temperature).

The sterilization instructions below have been validated to a sterility assurance level of 10⁻⁶.

1. Wrap instruments.
2. Inspect the packaging to ensure no rips, punctures, or seal failures are present in or on the packaging prior to loading into the sterilizer.
3. Load the instruments into the sterilizer by following the sterilizer manufacturer's recommended loading procedures and load configurations.
4. Follow the sterilizer manufacturer's recommended procedures to program the sterilizer with any one of the sets of sterilization cycle parameters found in Table 1 (for medical facilities in the US and its territories) and Table 2 (for medical facilities outside the US and its territories).

Table 1: Sterilization Cycle Parameters for the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	MINIMUM DRY TIME ¹
Steam	Gravity Displacement	250°F (121°C)	30 Minutes	30 Minutes
Steam	Gravity Displacement	270°F (132°C)	15 Minutes	30 Minutes
Steam	Gravity Displacement	275°F (135°C)	10 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	270°F (132°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	275°F (135°C)	3 Minutes	16 Minutes

For Medical Facilities Located Outside the United States and Its Territories: Some non-U.S. health care authorities recommend sterilization according to these parameters so as to minimize the potential risk of transmission of Creutzfeldt-Jakob disease, especially of surgical instruments that could come into contact with the central nervous system.

Table 2: Sterilization Cycle Parameters for Medical Facilities Outside the United States and Its Territories

METHOD	CYCLE	TEMPERATURE	EXPOSURE TIME	MINIMUM DRY TIME ¹
Steam	Gravity Displacement	273°F (134°C)	20 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	4 Minutes	30 Minutes
Steam	Dynamic-Air-Removal	273°F (134°C)	20 Minutes	30 Minutes

¹ The minimum dry times were validated using sterilizers having vacuum drying capabilities. Drying cycles using ambient atmospheric pressure may require longer dry times. Refer to the sterilizer manufacturer's recommendations.

NOTE: Chamber size and chamber load differences may exist between industrial and health care facility sterilizer models. The sterilization parameters listed in Tables 1 and 2 can be achieved in both health care facility and larger, industrial sterilizer models. Because of the many variables involved in sterilization, each medical facility should calibrate and verify the sterilization process (e.g. temperatures, times) used for their equipment.

The sterilization cycles listed in Table 2 are not considered by the FDA to be standard sterilization cycles. It is the user's responsibility to use only sterilizers and accessories (such as sterilization wraps, sterilization pouches, sterilization cassettes) cleared by the FDA for the selected sterilization cycle specifications (time and temperature).

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PRODUCT COMPLAINTS

Any health care professional (e.g., user of this system of products) having any complaints or having experienced any dissatisfaction in the product quality, identity, durability, reliability, safety, effectiveness, or performance should notify the distributor or Medtronic. Further, if any of the implanted spinal system component(s) ever "malfunctions" (i.e., does not meet any of its performance specifications or otherwise does not perform as intended), or is suspected of doing so, the distributor should be notified immediately. If any Medtronic product ever malfunctions and may have caused or contributed to the death or serious injury of a patient, the distributor should be notified immediately by telephone, fax, or written correspondence. When filing a complaint, provide the component(s) name and number, lot number(s), your name and address, the nature of the complaint, and notification of whether or not a written report from the distributor is requested.

© 2015 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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For a long time, epidemiologists have been aware that the epidemiology of hepatitis C virus (HCV) infection is different from that of hepatitis B virus (HBV) infection. The epidemiology of HCV infection is characterized by a high prevalence of chronic infection, a high rate of progression to cirrhosis and liver cancer, and a low rate of spontaneous clearance. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain populations, such as injection drug users, and a low prevalence of infection in other populations, such as health-care workers. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain countries, such as the United States, and a low prevalence of infection in other countries, such as Japan. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain regions, such as the Mediterranean basin, and a low prevalence of infection in other regions, such as the Amazon basin. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain ethnic groups, such as the Italian population, and a low prevalence of infection in other ethnic groups, such as the Chinese population. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain occupations, such as the military, and a low prevalence of infection in other occupations, such as the agricultural sector. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain social classes, such as the lower social class, and a low prevalence of infection in other social classes, such as the upper social class. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain countries, such as the United States, and a low prevalence of infection in other countries, such as Japan. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain regions, such as the Mediterranean basin, and a low prevalence of infection in other regions, such as the Amazon basin. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain ethnic groups, such as the Italian population, and a low prevalence of infection in other ethnic groups, such as the Chinese population. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain occupations, such as the military, and a low prevalence of infection in other occupations, such as the agricultural sector. The epidemiology of HCV infection is also characterized by a high prevalence of infection in certain social classes, such as the lower social class, and a low prevalence of infection in other social classes, such as the upper social class.

Abstract: The population of the city of Moscow is growing rapidly, and the number of people living in the city is increasing. The population of the city of Moscow is growing rapidly, and the number of people living in the city is increasing.

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NOT FOR DISTRIBUTION IN U.S.A. OR ITS TERRITORIES

Disassembly/Assembly Instructions for NFORCE™ Instruments

M708348B578 Rev. A

2016-07-22

Medtronic

* M 7 0 8 3 4 8 B 5 7 8 - A *

ENGLISH

IMPORTANT INFORMATION ON DISASSEMBLY/ ASSEMBLY FOR NFORCE™ INSTRUMENTS

Refer to the Medtronic Reusable Instruments Instructions for Use provided with the instruments for important information including reprocessing and sterilization instructions. Disassembly is required prior to reprocessing the instrument.

NFORCE™ Instruments are a set of instruments composed of distraction modules (Rack, Pinion, Key and Nut), lordosis modules (Rack and Ratchet), and extenders (Extender Open or Extender Adapter).

Description	Part Numbers
DISTRATOR RACK/DISTRATOR RACK L	5580101/5580111
DISTRATOR PINION/DISTRATOR PINION L	5580102/5580112
DISTRATOR KEY/DISTRATOR KEY L	5580103/5580113
DISTRATOR NUT/DISTRATOR NUT L	5580104/5580114
LORDOSIS RATCHET	5580115
LORDOSIS RACK	5580116
EXTENDER OPEN/EXTENDER ADAPTER	5580130/5580120

ASSEMBLY

The construction of NFORCE™ Instruments is completed in the Operating Room during the surgical procedure. The assembly instructions are included in the corresponding surgical technique document. See Figure 1 for a fully assembled NFORCE™ Instrument system.

Note: the system can be assembled with either the Extender Open or with Extender Adapter depending on surgery approach (Figure 1). The disassembly steps are the same in both cases. Figures 2 to 8 show only the Extender Open.

DISASSEMBLY

1. Remove Distractor Key (#5580103/#5580113). See Figure 2.
2. Unscrew Distractor Nuts (#5580104/#5580114) from Extenders (#5580130/#5580120). See Figure 3.
3. Unplug lordosis modules from Extenders (#5580130/#5580120). See Figure 4. The lordosis modules are the assembly of Lordosis Ratchet (#5580115) with Lordosis Rack (#5580116).
4. Unplug distraction module from Extenders (#5580130/#5580120). See Figure 5. The distraction module is the assembly of Distractor Pinion (#5580102/#5580112) with Distractor Rack (#5580101/#5580111).
5. Disassemble distraction module:
 - 5.1. Set the button to central position on Distractor Pinion (#5580102/#5580112). See Figure 6.
 - 5.2. Slide out Distractor Pinion (#5580102/#5580112) from Distractor Rack (#5580101/#5580111).

See Figure 6.

6. Disassemble lordosis modules:

6.1. Press the button of Lordosis Ratchet (#5580115). See Figure 7.

6.2. Slide out Lordosis Rack (#5580116) from Lordosis Ratchet (#5580115). See Figure 7.

Note: parts of NFORCE™ Instruments are placed unassembled into the spinal instrument tray to be sterilized. See Figure 8.

©2016 Medtronic Sofamor Danek USA, Inc. All rights reserved.

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2020-03-10

Medtronic

IMPORTANT INFORMATION ON THE TRANSLACE™ SPINAL TETHERING SYSTEM

PURPOSE

The Translace™ Spinal Tethering System is intended to be used in conjunction with pedicle screw and rod based posterior fixation systems that help provide immobilization and stabilization of spinal segments.

DESCRIPTION

The Translace™ Spinal Tethering System consists of a connector, a tether, and system specific instrumentation that allows the spine to be secured to a rod construct at levels where the use of a bone screw is not optimal. The Translace™ Spinal Tethering System acts as an alternative to sublaminar wires and hooks. The Translace™ Spinal Tethering System is intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures. The Translace™ Spinal Tethering System will supplement the interface between spinal anatomy and rod, through a pedicle sparing technique. The tether was designed to be used in conjunction with the connector. This system should be implanted only by a physician thoroughly knowledgeable in the implant's material and surgical aspects and instructed as to its mechanical and material applications and limitations.

INDICATIONS

The Translace™ Spinal Tethering System is a temporary implant for use in orthopedic surgery. The system is intended to provide temporary stabilization as a bone anchor during the development of solid bony fusion and aid in the repair of bone fractures. Indications for use include the following applications:

1. Spinal trauma surgery, used in sublaminar or facet wiring techniques.
2. Spinal reconstructive surgery, incorporated into constructs for the purpose of correction of spinal deformities such as idiopathic, congenital and neuromuscular scoliosis in patients 8 years of age and older, adult scoliosis, kyphosis, and spondylolisthesis.
3. Spinal degenerative surgery, as an adjunct to spinal fusions.

The Translace™ Spinal Tethering System may also be used in conjunction with other medical implants made of similar metals whenever "wiring" may help secure the attachment of other implants.

CONTRAINDICATIONS

- Bone metabolism disorders that potentially compromise the mechanical support expected for this type of implant (any abnormality affecting the normal functioning of bone tissue including acute spinal osteoporosis, bone resorption, osteopenia, primary or metastatic tumors of the spine, active infection at the site, and some metabolic disorders that affect osteogenesis).
- Insufficient quality and quantity of bone which would limit the efficacy of osteosynthesis.
- Severe fractures such that segments may not be maintained in satisfactory proximate reduction.
- Fractures that potentially compromise the mechanical support expected for this type of implant.
- Active local or systemic infections, or recent history of local or systemic infections that may jeopardize the outcome of the operation.
- Major local inflammation.
- Open wounds.
- Immunosuppressive diseases.
- Any other medical or surgical condition which would preclude the potential benefit of spinal implant surgery, such as the presence of a tumor, congenital abnormalities, elevation of sedimentation rate unexplained by other diseases, elevation of white blood count (WBC), or a marked left shift in the WBC differential count.
- Pregnancy.
- Sensitivity to implant materials or foreign bodies. If there is any suspected sensitivity to the materials used, the patient should have the appropriate tests before selection and implantation of the material.
- In any situation where implant use would interfere with anatomical structures or expected physiological performance, such as impinging on vital structures.
- Any other relative contraindication including:
 - Obesity. A person who is overweight or obese can overload the system leading to failure of fixation or breakage of the material.

- Excessive physical activity. Intense occupational level or activity level of the patient or a state of senility, mental illness, or other use of psychoactive substances. These conditions, along with others, may lead the patient to neglect surgeon recommendations, in turn leading to failure of the fixation or rupture of the material.
- Any neuromuscular deficit resulting in an unusual overload of the system during the consolidation period.
- Patients with insufficient muscle or tissue coverage of the operative site.
- Disease conditions shown to be safely and predictably managed without the use of fixation devices.

PRECAUTIONS

- The implantation of this type of implant should be performed only by experienced spinal surgeons with specific training in the use of this spinal system because this is a technically demanding procedure presenting a risk of serious injury to the patient.
- The information in these instructions is necessary but not sufficient for using this system. This information in no way takes the place of the professional judgment, expertise, and experience of the surgeon in patient selection, preoperative planning and the choice of implant size, knowledge of the anatomy and biomechanics of the spine, knowledge of the materials and understanding of the mechanical characteristics of the implants used, training and expertise in spinal orthopedic surgery, the use of accessory instruments for implantation, and patient commitment to follow an appropriate postoperative plan and having the expected postoperative exams.
- A successful result is not always achieved in every surgical case. This fact is especially true in spinal surgery where many extenuating circumstances may compromise results.
- The device system is not intended to be the sole means of spinal support. Use of this product without a bone graft or in cases that develop into a non-union will not be successful. No spinal implant can withstand body loads without the support of bone. In this event, bending, loosening, disassembly, and/or breakage of the device(s) will eventually occur.
- Before using the implant, it is absolutely essential to check the integrity of the packaging and to check the expiration date on the label, which guarantees that sterility has been maintained.
- Only physicians who are familiar with and trained on the techniques for using Tether Connector instruments are authorized to use them.
- Instruments must be checked before the procedure to be sure they are not worn or damaged.
- Before use, it is advisable to verify the instruments are intact and functioning completely properly.
- Surgeons should ensure they are not using instruments that could cause inappropriate tension on the spinal column or on the implants and must scrupulously follow the operative protocol described in the surgical technique brochure available from Medtronic. This means, for example, that surgeons must avoid injuring the patient from pressure exerted during in-situ repositioning of the instrument.
- To reduce the risk of breaking, implants should not be bent, folded, struck, or scratched with instruments unless the Translace™ Spinal Tethering System surgical technique specifies otherwise.
- Instruments should be used with extreme caution near vital organs, nerves, and blood vessels.
- Instruments can be reused after being reprocessed unless specified.
- When hypersensitivity is known or suspected, it is recommended to check the skin tolerance of the implant's materials before implantation.

WARNINGS

- Depending on the resistance of the patient's bone, do not apply an excessive force on the braid that could lead to fracture of the lamina and/or transverse process.
- Do not use if package is opened or damaged or if expiration date has passed.
- Do not use damaged implants.
- Never reuse an implant. Even if it seems to be intact, a previously used implant can have imperfections or defects that could reduce its lifetime.
- Resterilization of this implant is strictly prohibited. If a single-use product is reused, the performance, cleaning and sterility of the device are no longer assured. This can in particular result in failure of the procedure or risks of infection that can lead to death of the patient.
- It is essential to adhere to aseptic conditions when opening the protective packaging and extracting the implant.
- It is extremely important to handle implants carefully. The surgeon and the surgeon's assistants should avoid nicking or scratching components.
- All implants should be used in the original form unless specifically stated. If applicable, any modification of the implant is exclusively the surgeon's responsibility.
- Only proper use of the specific accessory equipment for the implant ensures satisfactory implant placement. Before use, verify instruments are intact and functioning completely properly.
- The Translace™ Spinal Tethering System tether and connector were designed to be used together and not separately.
- Translace™ Spinal Tethering System has not been tested for use with other manufacturers products.
- The malleable tip of the tether is made from commercially pure titanium, or medical grade stainless steel, and polyester. This portion of the tether serves only to put the implant in place and should be removed once the connector is positioned and is correctly placed under tension. This removal is performed by cutting the tip and extra tether and is described in the surgical technique. This portion of the tether is not intended to be implanted.

- Medical grade titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy may be used together. Never use titanium, titanium alloy, and/or medical grade cobalt-chromium-molybdenum alloy with stainless steel in the same construct.

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The implantation of tethering systems in pediatric patients should be performed only by experienced spinal surgeons with specific training in the use of this tethering system in pediatric patients because this is a technically demanding procedure presenting a risk of serious injury to the patient.

Preoperative and operating procedures, including knowledge of surgical techniques, good reduction, and proper selection and placement of the implants are important considerations in the successful use of the system in pediatric patients.

Physician note: although the physician is the learned intermediary between the company and the patient, the important medical information in this document should be conveyed to the patient.

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IMPLANT SELECTION

Surgical implants are subject to repeated stresses in use, and their strength is limited by the need to adapt the design to the size and shape of human bones. Unless great care is taken in patient selection, proper placement of the implant, and postoperative management to minimize stresses on the implant, such stresses may cause implant fatigue and consequent breakage, bending, or loosening of the device before the healing process is complete which may result in further injury or the need to remove the device prematurely.

DEVICE FIXATION

The Translace™ Spinal Tethering System contains tethers and connectors to be used for implantation with cleared rods and implants of various diameters. Rods ranging from 4.5mm to 6.30mm, including CD Horizon™ Apex Rods and rods made of titanium alloy or cobalt chrome may be used with the Translace™ Spinal Tethering System. The spinal tether consists of a biocompatible polyester tether which connects through a titanium alloy connector around the lamina that is then fixated to the rod. Once the tether connector is implanted and the rod is placed, the tether connector is attached to the rod, the tethers are sequentially tightened, and the set screws are broken off. For self-breaking set screws, always hold the assembly with the counter torque device. Tighten and break-off the head of the set screw to leave the assembly at optimum fixation security. After the upper part of the self-breaking set screw is sheared off, further re-tightening is not necessary and not recommended. The head part should not remain in the patient. After the upper part of the self breaking set screw is sheared off, re-adjustment is not possible unless the set screw is removed and replaced with a new one.

PREOPERATIVE

- Only patients that meet the criteria described in the indications should be selected.
- Patient conditions and pre dispositions such as those addressed in the contraindications should be avoided.
- Care should be used when handling and storing implant components. Implants should not be scratched or otherwise damaged. Implants and instruments should be protected during storage, especially from corrosive environments.
- An adequate inventory of implants should be available at the time of surgery; normally a quantity in excess of what is expected to be used.
- Since mechanical parts are involved, the surgeon should be familiar with the various components before using the equipment and should personally assemble the devices to verify all parts and necessary instruments are present before surgery. Translace™ Spinal Tethering System components (described in the DESCRIPTION section) are not to be combined with components from another manufacturer.
- Instruments should be cleaned and sterilized before use. Additional sterile implants and instruments should be available in case of an unexpected need.

INTRAOPERATIVE

- Extreme caution should be used around the spinal cord and nerve roots. Damage to nerves will cause loss of neurological functions.
- Breakage, slippage, or misuse of instruments or implant components may cause injury to the patient or operative personnel.
- Metallic tether connectors should not be bent. Use great care to ensure implant surfaces are not scratched or notched since such actions may reduce the functional strength of the construct.
- Use an imaging system to facilitate surgery.
- Excessive force on the polyester tether during implantation and corrective procedures could lead to fracture of the lamina and/or transverse process.
- After set screw break-off, and before wound closure, the excess tether should be cut from the connectors, leaving a 5mm minimum tail.
- Bone graft must be placed in the area to be fused and graft material must extend from the upper vertebrae to the lower vertebrae being fused.
- Bone cement should not be used because safety and effectiveness of bone cement has not been determined for spinal uses, and this material will make removal of the components difficult or impossible. Heat generated from the curing process may also cause neurologic damage and bone necrosis.

- Before closing the soft tissues, provisionally tighten the nuts or screws, especially screws or nuts that have a break-off feature. Once this is completed, go back and firmly tighten all screws and nuts. Recheck the tightness of all nuts or screws after finishing to ensure none loosened during the tightening of the other nuts or screws. Failure to do so may cause loosening of the other components.

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- Detailed instructions on the use and limitations of the device should be given to the patient. If partial weight-bearing is recommended or required prior to firm bony union, the patient must be warned that bending, loosening, and/or breakage of the device(s) are complications which may occur as a result of excessive or early weight-bearing or muscular activity. Risk of bending, loosening, or breakage of a temporary internal fixation device during postoperative rehabilitation may be increased if the patient is active, debilitated, or demented. The patient should be warned to avoid falls or sudden jolts in spinal position.
- To allow the maximum chances for a successful surgical result, the patient or devices should not be exposed to mechanical vibrations or shock that may loosen the device construct. The patient should be warned of this possibility and instructed to limit and restrict physical activities, especially lifting and twisting motions and any type of sport participation. The patient should be advised not to smoke tobacco, use nicotine products, or consume alcohol or non-steroidals or anti-inflammatory medications such as aspirin during the bone graft healing process.
- The patient should be advised of their inability to bend or rotate at the point of spinal fusion and taught to compensate for this permanent physical restriction in body motion.
- Failure to immobilize a delayed or non-union of bone will result in excessive and repeated stresses on the implant. By the mechanism of fatigue, these stresses can cause the eventual bending, loosening, or breakage of the device(s). It is important that immobilization of the spinal surgical site be maintained until firm bony union is established and confirmed by roentgenographic examination. If a state of non-union persists or if the components loosen, bend, or break, the device(s) should be revised or removed immediately before serious injury occurs. The patient must be adequately warned of these hazards and closely supervised to ensure cooperation until bony union is confirmed.
- As a precaution, before patients with implants receive any subsequent surgery (such as dental procedures), prophylactic antibiotics may be considered, especially for high-risk patients.
- The Translace™ Spinal Tethering System implants are temporary internal fixation devices. Internal fixation devices are designed to stabilize the operative site during the normal healing process. After the spine is fused, these devices serve no functional purpose and may be removed. While the final decision on implant removal is, of course, up to the surgeon and patient, in most patients, removal is indicated because implants are not intended to transfer or support forces developed during normal activities. If the device is not removed following completion of its intended use, one or more of the following complications may occur: (1) corrosion, with localized tissue reaction or pain; (2) migration of implant position possibly resulting in injury; (3) risk of additional injury from postoperative trauma; (4) bending, loosening, and breakage which could make removal impractical or difficult; (5) pain, discomfort, or abnormal sensations due to the presence of the device; (6) increased risk of infection; (7) bone loss due to stress shielding; and (8) potential unknown or unexpected long term effects such as carcinogenesis. Implant removal should be followed by adequate postoperative management to avoid fracture, re-fracture, or other complications.
- Retrieved devices should be treated in such a manner that reuse in another surgical procedure is not possible. As with all orthopedic implants, the Translace™ Spinal Tethering System implants should never be reused under any circumstances.

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PACKAGING

If devices are individually packaged, the packages for each of the implants and/or instruments should be intact upon receipt. If a loaner or consignment system is used, all sets should be carefully checked for completeness and all components including instruments should be carefully checked to ensure there is no damage prior to use. Damaged packages or products should not be used and should be returned to Medtronic.

STERILE IMPLANTS

Implants are provided sterile and should only be used if they are marked sterile and clearly labeled as such in an unopened sterile package provided by Medtronic. Only sterile products should be placed in the operative field. Implants should never be reprocessed.

PRODUCT COMPLAINTS

To report any product problems, contact Medtronic.

MRI INFORMATION

The Translace™ Spinal Tethering System has not been evaluated for safety and compatibility in the MR environment. It has not been tested for heating, migration, or image artifact in the MR environment. The safety of the Translace™ Spinal Tethering System in the MR environment is unknown. Scanning a patient who has this device may result in patient injury.

FURTHER INFORMATION

Recommended directions for use of this system (surgical operative techniques) are available at no charge upon request. If further information is required, contact Medtronic.

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EXPLANATION OF SYMBOLS

CAUTION: Federal law (USA) restricts these devices to sale by or on the order of a physician.

Do not re-use

Batch code

Manufacturer

Catalogue number

For US audiences only

Sterilized using Irradiation

Use-by date

Consult instructions for use at this website.

Medical device

Single sterile barrier system

Double sterile barrier system

Single sterile barrier system with protective packaging outside

Single sterile barrier system with protective packaging inside

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Never use instruments or Medtronic Transportation/Sterilization Cassettes with obvious signs of excessive wear, damage, or that are incomplete or otherwise non-functional.

Medtronic Transportation/Sterilization Cassettes must be examined by the user prior to use. Examination should be thorough and must include a visual and functional inspection. Examination must also ensure the Medtronic Transportation/Sterilization Cassettes are complete. Do not proceed with reprocessing of a damaged Medtronic Transportation/Sterilization Cassette.

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этапами обработки не должно быть задержек. /Cleaning is the initial step. Sterilization occurs later in reprocessing and is intended to kill microorganisms to reduce the likelihood of transmission and possibility of infection. To ensure acceptable reprocessing, there should be no delay between the steps in this document.

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Instruments may be safely and effectively processed using manual or a combination of manual and automated washer- disinfecter (enzymatic or alkaline) cleaning procedures outlined in this document. Medtronic Transportation/Sterilization Cassettes may be safely and effectively processed using manual or a combination of manual and automated washer- disinfecter (enzymatic) cleaning procedures outlined in this document.

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Universal precautions for handling this device after use should be observed by all hospital personnel according to OSHA Standard 29 CFR 1910.1030, Occupational Exposure to Bloodborne Pathogens.

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Note: alkaline cleaning agents may require neutralization following cleaning. Refer to manufacturer's instructions to determine if neutralization is required, and to follow manufacturer's instructions for neutralization. It is critical that alkaline cleaning solutions be properly neutralized and rinsed from instruments.

Aluminum Transportation/Sterilization Cassettes: the use of neutral detergents, enzymatic cleaners (7.0 – 8.0 pH) soft bristled brushes and soft pipe cleaners are recommended. Do not use alkaline cleaning agents as the alkaline cleaners will cause significant damage and render to the aluminum Transportation/Sterilization Cassettes unusable.

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©2017 Medtronic Sofamor Danek USA, Inc. All rights reserved.
© ٢٠١٣ Gauthier Biomedical, Inc. All rights reserved
©Copyright 2013 Tecomet