SENSITEC

Negative Pressure Wound Therapy Devices SENSITEC With Accessories

EN

by Apexmed

**Negative pressure wound therapy devices SENSITEC,
available modifications: SensiVac WT-100, SensiVac WT-110,
SensiVac WT-120, SensiVac WT-200, SensiVac WT-300**

Maintenance Documentation

SENSITEC
by Apexmed

Contents

1	Introduction	5
2	Warnings	5
3	Safety Instructions	5
4	Clinical Warnings	6
5	Precautions	7
6	Specifications	9
7	SensiVac WT-100, SensiVac WT-110, SensiVac WT-120	12
7.1	Operation Instructions	14
7.2	Device Maintenance	16
7.3	Therapy Mode Choice and Settings	16
8	SensiVac WT-200, SensiVac WT-300	20
8.1	Operation Instructions	23
8.2	Device Maintenance	24
8.3	Therapy Mode Choice and Settings	24
9	Malfunctions and Troubleshooting	29
10	Service	30
11	Transportation and storage	30
12	Disposal	31
13	Contents	31
14	EMC Declaration	32
15	Limited Warranty	34
16	Symbols	35

1. Introduction

Apexmed International B.V. thanks you for the product purchase and offer to you to study the present user manual. The manual contains all necessary data on the device application. Follow the given instructions and keep the manual unchanged in a place accessible for operator. The manual's contents can be changed without preliminary notification and extra obligations, therefore the manufacturer is entitled to make changes and improvements. Copying and translation of the present manual without written permission of Apexmed International B.V. is forbidden.

Attention! Study the present manual before the device application.

Attention! The device is intended for application only on doctor's orders and under medical personnel control.

Attention! Never open the device, appeal to Apexmed International B.V. technical department if any problem occurs.

Attention! Apexmed International B.V. is not liable for any side effects resulted from erroneous actions considering device application and its accessories, recommended by the company.

2. Warnings

Read and study the present manual carefully before the device application. The device should be maintained according to the manual. Apexmed International B.V. is not liable for any incidental, collateral or indirect damages, incurred by the customer as a result of wrong application, improper maintenance or unauthorized repair.

3. Safety instructions

For correct device operation the following conditions are recommended. Disregard of the conditions leads to the set warranty termination.

- Use the device only according to the present user manual and labeling.
- Assure that main power of the room meets applicable national standards.
- In spite of the fact that the negative pressure wound therapy devices SENSITEC (NPW device SENSITEC) meet requirements on electromagnetic compatibility, sometimes electrical disturbances may occur in the devices. If there are electrical disturbances, replace the device in another place further from electrical emission source.
- For the device operation and accumulator charge use only power unit supplied by Apexmed International B.V. Usage of another power unit may lead to the device damage.
- Power cords may present a tripping hazard. Ensure that all cords are out of areas where people may walk.
- Do not operate this product if it has a damaged power cord, power supply or plug. If these components are worn or damaged, contact Apexmed International B.V. or an authorized provider.
- Do not drop or insert any object into any port or tubing of this product.
- Do not connect this product or its components to devices not recommended by Apexmed International B.V.
- Keep this product away from heated surfaces.
- Do not reach for a product that has fallen into water. Unplug the unit immediately if plugged into electrical source. Disconnect the unit from dressing and contact Apexmed International B.V. or an authorized provider.

Warning/Caution symbols used in this manual apply to hazardous or unsafe practices which could result in personal injury or property damage.

Attention! It is dangerous for patients to analyze the wound on their own according to therapy outcome. Clinician's advice should be strictly followed.

Attention! The device is not for use by infants or by those unable to express their thoughts. The device should not be used for any other purpose other than the intended use. All disassembly, adjustment, modification and/or repair are not permitted.

Warning! Do not strike or drop the device. Keep top side up when using. The device contains laser emission source. Do not watch into the holes on the right panel to prevent eyes injury.

The device is used in combination with wound dressing kits and a disposable canister. Ensure that the wound bed, the canister and other tubing are well sealed when using, or it may affect clinical performance of the device, and even cause device failure.

4. Clinical Application

Negative Wound Pressure device SENSITEC (Device) is an AC-powered portable device with battery back-up, contains portable vacuum pump. The device is used with disposable wound dressing kit, canister set to apply localized vacuum to draw wound edges together while removing exudate and promoting microcirculation and granulation tissue formation, microbial contamination reduction and protection from cross-infection. Exudate is collected from wound into disposable canister.


Fig. 1.

Wound treatment with negative pressure (vacuum therapy)

Healing of wounds by negative pressure – represents the innovative method facilitating the fast wound healing by building up continuous and/or intermittent negative pressure in the wound bed. Therefore the local interstitial edema of tissues is reduced, the local lymph circulation and transcapillary transfer is increased, which resulted in improvement of traumatic medium, lymphotrophy, and increasing the rate of granulations formation, and improvement of the wound bed perfusion additionally helps in the wound decontamination. Direct impact of the negative pressure on the wound bottom and edges under conditions of the external insulation exercises the permanent effect on the wound edges helping in the wound edges approximation. Such effect directly reduces the wound sizes, regardless the intensity of the cellular proliferation.

The abovementioned therapy mode has many indications, and very effective for healing of long lasting open wounds, as well as acute wounds.

The user manual on NPWT device SENSITEC is not intended as a guarantee of positive results of wound treatment. The manual is a guidance on the choice of the most efficient wound treatment protocol for each patient.

1. NPWT device SENSITEC is a hospital care device.
2. NPWT device SENSITEC must be operated by qualified physicians, nurses or healthcare personnel who have received training related to NPWT. Patients are forbidden to operate the device.
3. Accessories used with NPWT device SENSITEC must be certified and registered in Russian Federation.
4. NPWT device SENSITEC is not MRI, X-ray, CT, fluoroscopy or ultrasound compatible.
5. NPWT device SENSITEC is a medical device that requires special safety precautions and must be installed and placed in operation in accordance with the attached EMC information. Portable and mobile RF communication devices (mobile telephones) can affect the device operation.
6. NPWT device SENSITEC is not suitable for use in the presence of a flammable anaesthetic mixture with air or with oxygen or nitrous oxide.
7. NPWT device SENSITEC is not suitable for use in areas where there is danger of explosion (e.g. hyperbaric oxygen unit).
8. Patients must be closely monitored for bleeding. If sudden or increased bleeding is observed, immediately discontinue therapy, take appropriate measures to stop bleeding.
9. Patients suffering from difficult hemostasis or who are receiving anticoagulant therapy have an increased risk of bleeding. During therapy, avoid using hemostatic products that may increase the risk of bleeding, if disrupted.
10. Sharp edges or bone fragments in a wound must be covered or removed prior to using NPWT device SENSITEC due to risk of puncturing organs or blood vessels while under negative pressure.
11. Do not use NPWT device SENSITEC on exposed blood vessels or organs.
12. In the event defibrillation is required, disconnect the device from the wound dressing prior to defibrillation. Remove the wound dressing if it will interfere with defibrillation.
13. When operating, transporting, repairing or disposing of NPWT device SENSITEC and its accessories, the risk of infectious liquids being aspirated or contamination of the device through incorrect use cannot be eliminated. Universal precaution should be observed whenever working with potentially contaminated parts or equipment. To avoid contact with blood and other biological liquids personal protection equipment must be used (disposable medical gloves).
14. Do not use on infants or patients who are unable to express their thoughts.
15. Do not assemble, adjust, modify or repair by yourself. Please refer to Apexmed International B.V. authorized agent for service.
16. Wound dressing kit, canister, tubing, and bacteria filter are disposable, and must not be reused.

5. Precautions

1. Precautions should be taken for patients who are or may be:
 - Receiving anticoagulant therapy or platelet aggregation inhibitors, actively
 - bleeding or have weakened blood vessels or organs.
 - Suffering from difficult wound hemostasis.
 - Untreated for malnutrition.
 - Noncompliant or combative.
 - Suffering from wounds in close proximity to blood vessels or delicate fascia.
2. Infected wounds may require more frequent dressing changes. Regular monitoring of the wound must be maintained to check for signs of infection.
3. Duration of vacuum therapy cycle is determined for each patient individually, taking into consideration wound location and state, exudate quantity, dressing seal integrity, wound contamination, etc.
4. This device should only be used by qualified and authorized personnel obtaining substantial knowledge on vacuum wound therapy.
5. If the device has been at temperatures below freezing, the device must be brought to room temperature for at least 12 hours prior to use or the pumping unit may be damaged.
6. Ensure the tubing is installed completely and without any kinks to avoid leaks or blockages in the vacuum circuit.
7. Position the device and tubing appropriately to avoid the risk of causing a trip hazard and so the patient does not lie on the tubing.
8. When bathing/showering, the patient must be disconnected from the device.
9. If any liquids penetrate the pump, discontinue use and return to Apexmed International B.V. authorized provider for service.
10. Check the in-line overflow protection and replace as necessary.

11. Underlying structures, such as bone and tendon, must be covered by the no adherent dressing layer.
12. Inspect the dressing site as appropriate based on therapeutic need.
13. The vacuum level should never be painful. If the patient reports discomfort with the vacuum level, consider reducing the pressure.
14. Ensure dressing clamp is engaged before switching off the device. For reconnection to the device, ensure that the device is active prior to opening the dressing clamp.

Indications for Use

Negative Pressure Wound Therapy Device SENSITEC is applied with standard wound dressing kits for various acute and chronic wounds treatment.

- Burns of I and II degree;
- Surgical wounds;
- Trauma-induced wounds;
- Diabetic / neuropathic ulcers;
- Bedsores;
- Venous and arterial ulcers;
- Radiation ulcers.

Contraindications

- Necrotic tissue with eschar;
- Blood coagulation abnormalities with hemorrhage risk;
- Untreated osteomyelitis (relative);
- Malignancy in the wound (excluding palliative therapy for life quality improvement);
- Exposed vasculature, nerves, anastomotic site or organs;
- Non-enteric and unexplored fistulas.

6. Specifications

Main dimensions and parameters

Model	SensiVac WT-100
External view	
Vacuum range	- 55 MmHg ~ - 225 MmHg
Therapy mode	Continuous, intermittent
Max flow	3 l/min.
Battery Run Life	average 40 hours (under the vacuum -155 MmHg, at least)
Laser module for overflow protection	No
Canister	800 ml
Dimension	270x228x205 mm
Weight	2,4 kg
Model	SensiVac WT-110
External view	
Vacuum range	- 55 MmHg ~ - 225 MmHg
Therapy mode	Continuous, intermittent
Max flow	3 l/min.
Battery Run Life	average 20 hours (under the vacuum -155 MmHg, at least)
Characteristics of laser module for overflow protection	
Light source	LED
Emission wave length	660 nm
Laser emission maximal output power	1 mW max
Canister	800 ml
Dimension	270x228x205 mm
Weight	2,4 kg

Model	SensiVac WT-120
External view	
Vacuum range	- 55 MmHg ~ - 225 MmHg
Therapy mode	Continuous, intermittent, combined
Max flow	3 l/min.
Battery Run Life	average 20 hours (under the vacuum -155 MmHg, at least)
Characteristics of laser module for overflow protection	
Light source	LED
Emission wave length	660 nm
Laser emission maximal output power	1 mW max
Canister	800 ml
Dimension	270x228x205 mm
Weight	2,4 kg
Model	SensiVac WT-200
External view	
Vacuum range	- 55 MmHg ~ - 225 MmHg
Therapy mode	Continuous, intermittent
Max flow	4 l/min.
Battery Run Life	average 10 hours (under the vacuum -155 MmHg, at least)
Sensor display	
Characteristics of laser module for overflow protection	
Light source	LED
Emission wave length	660 nm
Laser emission maximal output power	1 mW max
Canister	350 ml
Dimension	265x200x102 mm
Weight	2,0 kg

Model	SensiVac WT-300
External view	
Vacuum range	- 55 MmHg ~ - 225 MmHg
Therapy mode	Continuous, intermittent
Max flow	4 l/min.
Battery Run Life	average 10 hours (under the vacuum -155 MmHg, at least)
Sensor display	
Characteristics of laser module for overflow protection	
Light source	LED
Emission wave length	660 nm
Laser emission maximal output power	1 mW max
Canister	350 ml
Dimension	265x200x102 ml
Weight	2,0 kg

Time of operation mode setting – 20 sec. max

Noise level 60 dB max.

Negative Pressure Wound Therapy Device SENSITEC, available modifications: SensiVac WT-100, SensiVac WT-110, SensiVac WT-120, SensiVac WT-200, SensiVac WT-300 must be used in dry conditions only, protected from atmospheric effects.

- Working temperature: from +10 °C to + 40°C
- Relative air humidity: from 20 % to 85%
- Atmosphere pressure 700 gPa to 1060 gPa

Electrical characteristics

Power	Accumulator battery; AC current, with AC current adapter
Rated voltage	18 V DC
Maximum operating current	1,4 A
Adapter characteristics:	
Power source	External power unit (adapter)
External power unit voltage	(100 – 240) V
Nominal frequency	47/63 Hz
Input current	(1,62 – 0,72) A
Grounding	No
Output power unit voltage	18 V DC
Maximal output current strength of power unit	3,5 A
Battery characteristics	
Type	Lithium polymer battery
Maximum operating current	11,1 V DC
Capacity	3800 mAh
Charge time	4 h
Battery run life	40 h max

Reliability requirements

1.	Depending on potential impact of failure in the course of usage	Class B
2.	Water protection degree	IPX1
3.	Stability depending on sustained mechanical actions	For group 2

Settings

Settings	Time, therapy mode, functions
Therapy modes	Intermittent
	Continuous
	Combined (intermittent+continuous)
Exudate control	Exudate volume control
Alarms and warnings	Visual and auditory alarms
Data record display	Auto lock out function
	Total record of device operation data
Monitor and Indications	Key components monitoring & protection
	Pressure monitoring & protection
	Battery monitoring & protection
	Working environment temperature monitoring

All accessories are intended for single use only!

The aspirator for vacuum wound therapy Sensitec is equipped with integrated software version: 2.20140715.2, date of issue - 15/07/2014. The software is activated automatically after the device is switched on. It is used for the following functions:

- Selection and adjustment of therapy mode
- Therapy process visualization on a display in real time mode

The software structure is represented by the menu with the following parts:

- Self-testing
- Setup data input
- Archive (internal archive, statistics)

Protection level of the software from purposive and unpurposed changes corresponds to C group. Access to functions of software setting changes is limited for common users.

7. SensiVac WT-100, SensiVac WT-110, SensiVac WT-120

Panel & Keypad

Fig. 2. Device front panel и main menu interface.

EXPLANATION OF BUTTONS

Button	Function	Remarks
	Power on/off	Press down the button for at least one second.
Menu	Main menu	You cannot return to the main menu when therapy is on.
Lock	Key Lock/Unlock	Only when therapy is on can the keypad be locked. Press the LOCK button to lock keypad; Unlock by pressing the MENU button first, then the LOCK button.
Mute	Mute on/off	The default setting is mute off.
F1 / F2 / F3	Select/Enter/Skip	The function of the buttons changes according to the specific page.
	Select/Enter/Skip	The function of the buttons changes according to the specific page.
	Select/Enter/Skip	The function of the buttons changes according to the specific page.
	OK/Enter/Skip	The function of the buttons changes according to the specific page.

Warning! When therapy is on, keeping the device muted may cause delayed handling of risks for caregivers.

LED Description


ON Left	Power Indicator	The power indicator light is constantly illuminated while the device is on. The power indicator light is off while the device is off.
ON Right	Alarm & warning indicator	Refer to Troubleshooting for details.

Display Panel Description

Therapy ON	Functional display area for key F1. Press F1, The device will operate in the mode showed on display. If this display area is clear it shows that F1 function is not available. If display shows "Therapy ON" treatment starts after F1 button is pressed.
Intermittent	Functional display area for key F1. Press F1, The device will operate in the mode showed on display. If this display area is clear it shows that F1 function is not available. If display shows "Therapy ON" treatment starts after F1 button is pressed.
Continuous	Functional display area for key F3. Press F3, The device will operate in the displayed mode. If this display area is clear it shows that F3 function is not available. If display shows "Continuous" pressing of F3 button will make available interface of continuous mode settings.
Help	Functional display area for key  Press  The device will operate in the displayed mode. If this display area is clear it shows that function is not available. If display shows "Help" pressing of  button will make available interface of help mode settings.


7.1. Operation Instructions/User Manual

Attach the power

Plug the AC cord plug into an AC power adapter.

Fig. 3

Plug the DC connector into the power connection on the back of the device.

Fig. 4

Warning! External power unit is an integrated part of the device. Use only the power unit, supplied with NPWT SENSITEC device. Usage of another power unit may cause device damage.

Warning! Connect the device to mains power with plug. Main power source should be easily accessible to avoid problems during device operation.

Warning! Outlet ends of the power unit should be inaccessible for patient. Operator should avoid simultaneous touch of patient and outlet ends of the power unit.

Negative pressure wound therapy devices SENSITEC, models: SensiVac WT-100, SensiVac WT-110, SensiVac WT-120

Device preparation to work

NPWT SENSITEC device should be used with canister and wound dressing kit recommended by Apexmed International B.V.


Fig. 5

- 1 Device
- 2 Vacuum aspiration port
- 3 Connection tubing
- 4 Canister vacuum port
- 5 Container drainage port
- 6 Drainage tube
- 7 Canister
- 8 Wound dressing

Fig. 6

Place the canister holder into the slot on the right of the unit.

Ensure that there is no clearance between the canister holder and the device.


Fig. 7

Close the canister lid firmly, be assured that the lid clamps (6 clamps) are situated behind the canister edge.


Fig. 8

Put the canister vertically into the canister holder and twist to ensure that it locks in place. Adjust the position of the canister. Make the vacuum port close to pump-canister connector and the patient port outward.


Fig. 9

Connect one end of the tubing to aspirator vacuum port.


Fig. 10

Connect the other end of the tubing to canister vacuum port.


Fig. 11

Connect a separate connector from the kit to patient's port on the canister.


Fig. 12

After proper dressing kit application on the wound connect drainage tube from the kit to the canister drainage port (see User Manual on Dressing Kit Application). Tube blockage (e.g., tube kinking or sediment blockage) should be avoided. Aspirator is ready to perform therapy in a set mode (see Section Device Application).

Canister and Connection Tubing Replacement

1. Assure the vacuum is turned off.
2. Remove the drainage tubing between the canister and wound dressing, taking care to avoid the chance of fluid escaping. Close patient port with attached caps.
3. Disconnect the connective tubing between the canister vacuum port and the device vacuum port.
4. Take the canister out of the holder.
5. Install a new canister as described in instruction.

Attention! Integral parts of canister set are for single use only. One kit can be use for one patient. Before canister replacement assure that the canister, its lid or connection tubing are not damaged. Otherwise they should be replaced. It is recommended to replace canister once a week. Application of antibacterial filter between device and canister is recommended.

Remove or Replacement of Wound Dressing

1. Turn the vacuum off and wait until the vacuum pressure drops to 0 mmHg.
2. Disconnect the drainage tube from the canister, remove film cover, fixing drainage tube port, and remove port together with tube from the wound with foam or gauze bandage, trying to avoid exudate leakage.
3. Remove bandage out of wound.
4. Apply new bandage from wound dressing kit following user instruction.

Attention! If bandage is sealed and patient feels comfortable, the bandage replacement is held 2-3 times a week. If there is too much exudate more frequent replacement can be required. Infected wounds also require more frequent bandage replacement. If there are signs of progressive inflammatory process, appeal to doctor immediately.

Attention! Disposal of used canisters and tubing must follow facility protocols relating to the handling of potentially infected or biohazardous materials. Follow local governing ordinances regarding disposal of canisters and tubing.

7.2. Device Maintenance

Power ON/OFF

Fig. 13. Main screen view.

Power On: Continue to press power button until the following screen appears (see Fig. 13).

Power Off: continue to press powerbutton  for more than 1 sec. and release when the screen is off.

Attention! Do not force shutdown when therapy is on. Do not shut off the device until the vacuum pressure drops to 0.

7.3. Therapy Settings

Available operation modes for the device are: continuous and intermittent. User can choose the required mode by himself. Settings in the shadowed areas can be adjusted by buttons ▲ ▼.

Intermittent Therapy


Fig. 14.

1. On the main menu screen, press F2 to select <Intermittent> and enter into intermittent therapy setting (Fig. 14).
2. Press F1 to set therapy time. Use ▲ and ▼ to reduce or increase required time.
3. Press F2 to set vacuum pressure level. Use ▲ and ▼ to increase or reduce required vacuum pressure.
4. Press OK to confirm settings and access to main menu.
5. During the intermittent therapy setting, pressing F3 at any time can initiate the therapy with the current therapy settings.

Attention! On the intermittent therapy setting screen pressing the MENU button before pressing OK or F3 leads to making the changes invalid.

Treatment time could be set for 1-12, 24, 48, 72 hours, device stops when treatment is finished. Setting the treatment time at ∞ hour will make the treatment time indefinite. Therapy will not stop unless user press F1 (← pause). The therapy remaining time LED shows N/A and the operating time will indicate how long the device works.

Progressive Settings of Intermittent Therapy Mode

Fig. 15.

Fig. 16.

1. In the main menu after choosing [Help] – [System] – [Progressive] – [Intermittent], progressive settings screen will have the following view (Fig. 15).
 2. If [Inactive] is chosen, progressive settings of intermittent therapy mode remain on default (negative pressure impact time is 5 min, vacuum absence time – 2 min.)
 3. If vacuum impact interval should be changed, choose [Active] and enter required settings of intermittent therapy mode – negative pressure impact time and vacuum absence interval within one cycle (Fig. 16).
 4. Progressive settings of intermittent therapy mode (Fig. 16).
 5. Press F1 to determine negative pressure impact interval. Use ▲ and ▼ to reduce and increase interval.
- Press F3 to choose vacuum absence interval. Use ▲ and ▼ to increase and reduce interval.
- Press OK, to confirm the settings and return to main menu screen.

Continuous Therapy Mode


Fig. 17.

1. Press F3 on the front panel of the device to choose <Continuous>, and enter continuous therapy mode settings (Fig. 17).
2. Press F1, to set therapy time. Use ▲ and ▼ to increase or reduce required time.
3. Press F2 to set vacuum pressure level.
Use ▲ and ▼ to increase or reduce required vacuum pressure.
4. Press [OK] to confirm settings and enter the main menu.
5. Therapy with the chosen settings can be started after entering continuous therapy mode (F3) by pressing F3 button.

Attention! On the continuous therapy setting screen pressing the MENU button before pressing [OK] or F3 leads to making the changes invalid.

Treatment time could be set for 1-12, 24, 48, 72 hours, device stops when treatment is finished. Setting the treatment time at ∞ hour will make the treatment time indefinite. Therapy will not stop unless user press F1 (⏸- pause). The therapy remaining time LED shows N/A and the operating time will indicate how long the device works.

Combined therapy mode


Fig. 17a.

1. Press on the front panel of device F2 and F3 buttons at once to choose "Combined" and enter the therapy mode setting.
2. Set parameters of intermittent therapy (see above).
3. Choose time interval between therapy modes by pressing F1 button.
4. Set parameters of continuous therapy (see above).

Therapy On

1. Press <Menu> button on the front panel to enter the main menu.
2. Press F1 button <Therapy On>, to start therapy according to recently chosen settings (See Section Therapy Settings).
3. During intermittent/ continuous therapy setting, pressing F3 at any time can initiate the therapy.
4. If the progressive settings are active in intermittent therapy mode, active vacuum time can be set between 2 and 5 min. with 1 min. interval, and no vacuum time can be set between 1 and 3 min. with 1 min. interval. On default settings: Vacuum is on -5 min. No vacuum - 2 min.
5. Therapy must be stopped before skipping to the main menu.
6. Once the treatment time has ended, the system will automatically stop the therapy.
7. The keypad will be automatically locked if there is no operation within 1 minute. To unlock the keypad, press <MENU> first, and then press <LOCK>.
8. Treatment can be stopped / restarted anytime by pressing F1 button (⏸-pause/ II – start).
9. If it is necessary to change therapy mode or settings, return to the main menu and choose new settings as described in Section Therapy Settings.


Fig. 18.


Fig. 19.

The table contains explanations of word and graphic messages on the device display during intermittent (Fig. 18) and continuous (Fig.19) therapy.

No.	Comments
1	Displays therapy mode.
2	It displays Off when therapy is on, On when therapy pauses or stops.
3	Treatment progress bar. It indicates therapy progress.
4	Sealing degree. The shading indicates the sealing degree. More shading means that the system leakage is higher.
5	The set value of negative pressure (mm Hg)
6	Symbol of therapy mode (intermittent/continuous)
7	Current device state LED (ON/OFF)
8	The real-time vacuum pressure (mm Hg)
9	Exudate volume data*.

* - Canisters should meet requirements for the device.

Fig. 20.

Fig. 21.

Information Review

Contents	Operating Instructions
Time and date and their setting	Choose [Help]-[Time] in the main menu with corresponding buttons (F1, F2, <▲>, <▼>) set required values. To return to the menu [Help] press <  >.
Information on power supply source, component temperature, aspirator model, progressive settings of intermittent therapy	Operating Instructions <  >.
Contact data	Choose [Help]-[Contact] in the main menu.
Recent therapy cycle data	Choose [Help]-[Contact] in the main menu. To return to the menu [Help] press <  >.
Exudate volume data	Choose [Help]- [Exudate] in the main menu. To return to the menu [Help] press <  >.

8. SensiVac WT-200, SensiVac WT-300

Device Components and Control Panel


Fig. 22. Front view.

- 1 Sensor display
- 2 Power ON/OFF button
- 3 Button of canister detachment
- 4 Canister
- 5 LED


Fig. 23. Rear view.

- 1 Hidden handle
- 2 Holders

Device interface in "Doctor" mode


Fig. 24. Main menu view in "Doctor" mode

- 1 Warnings display area
- 2 System time
- 3 Sealing area
- 4 Current vacuum value
- 5 Exudate volume
- 6 Device operation state
- 7 Sensor buttons (icons)
- 8 Therapy progress
- 9 Prompt display area
- 10 Battery charge level

	Current mode is intended for doctor. To shift tot patient mode press icon
	System information screen access
	Press to check device records
	Help
	Therapy ON
	Therapy OFF (or stop – pause)
	Symbol blinking means that the therapy is over, symbol rotating means that the therapy is in progress
	Bad sealing

Note 1: Press icon "Exudate volume" to enter interface of exudate volume record.

Note 2: "Sealing level display area". If there is sealing degradation the icon displays more bubbles.

Device interface in "patient" mode


Fig.25. Main screen view in "patient" mode.

	Current mode is intended for patient. To enter "Doctor" mode press the icon and fulfill actions displayed on the screen.
---	--

8.1. Device Operation


Fig. 26. Device control is performed by sensor display.

Power ON

Insert pin plug into AC-adaptor. Connect power unit to side-panel of the device.

Device Preparation to Work

NPWT SENSITEC device should be used with canister and wound dressing kit recommended by Apexmed International B.V.

Step 1: Insert canister into device.

Step 2: Connect one end of the tubing to aspirator vacuum port.

Step 3: Connect the other end of the tubing to canister vacuum port.

Step 4: After proper dressing kit application on the wound connect drainage tube from the kit to the canister drainage port (see User Manual on Dressing Kit Application). Tube blockage (e.g., tube kinking or sediment blockage) should be avoided.

Canister and Connection Tubing Replacement

1. Switch off the device.
2. Lock all clamps of drainage tube as shown on figure.
3. Disconnect canister tubing from dressing drainage tubes.
4. Press button . Canister will be detached, but remains in slot.
5. Remove canister out of the slot.
6. Canister disposal should be performed according to the hospital's protocol.

Attention! If bandage is sealed and patient feels comfortable, the bandage replacement is held 2-3 times a week. If exudate is too dense or contains clots, more frequent replacement can be required. Infected wounds also require more frequent bandage replacement. If there are signs of progressive inflammatory process appeal to doctor immediately.

Attention! Used-up canisters and dressing materials disposal should correspond to medical institution's rules on potentially infected and biologically hazardous materials treatment. Local laws on canisters and dressing materials disposal must be followed.

8.2. Device Maintenance

Power ON

- Connect the device to mains power.
- Press and hold for 2 sec. button  to start the device .

Power OFF

Press and hold over 1 sec. button  to switch the device off.

8.3. Choice and setting of therapy mode

There are two modes available for the device: intermittent and continuous.


Fig. 27. Intermittent mode setting I.

Parameter setting	Note
	Intermittent mode time can be set between 1~12\24\36\72\∞ hours, to choose required interval use buttons  or  .
	Intermittent mode at low pressure can vary 1 min\ 2 min\ 3 min, to choose required value use buttons  or  .


Fig. 28. Intermittent mode setting II.

Parameter setting	Note
	Intermittent mode period can be set with the following interval 1~12\24\36\72\∞ hours, to choose required value use buttons ◀ or ▶ ∞ means no time limit, the device will work until user turn off the device, until battery charge is over, until automatic turn off of the device.
	Intermittent mode therapy with high pressure can be set between the following pressure values 55 mm Hg.~225 mm Hg, to choose required value use ◀ or ▶. Interval therapy at low pressure is constantly set on 0 mm Hg.

Description of Sensor Buttons (Icons)

	Shift to window "Intermittent therapy mode settings".
	Shift to window "Continuous therapy mode settings".
	Shift to window "Date and time setting".
	Return to "Doctor's home page" without setting save.
	Confirm and save all settings and return to "Doctor's home page".
	Press the button to start treatment.
	Press this button to shift to the next window.
	Press this button to return to the previous menu section.


Fig. 29. Continuous therapy mode setting.

Combined Therapy Mode Setting

Press and hold intermittent mode button and press at the same time continuous mode button to enter combined mode setting.

Initially after entering the mode set. With buttons or set required start/stop time of treatment (in minutes).

Use buttons ◀ or ▶ and to change pressure value. Set pressure range ◀ or ▶ lies within – 55 mm Hg... – 225 mm Hg.

Press this button to save intermittent therapy settings and entering continuous therapy parameters in combined therapy mode.

With buttons and required total therapy time can be set (in hours).

Use buttons and to change pressure.

Set pressure value range lies within – 55 mm Hg ... – 225 mm Hg with intervals 10 mm Hg.

Parameter setting	Note
 ∞ Hour 	Therapy time for continuous mode can vary within 1~12\24\36\72\∞ hours, to choose required value use  or  buttons. ∞ means no time limit, the device will work until user turn off the device, until battery charge is over, until automatic turn off of the device
 55 mmPn 	Pressure value during continuous therapy can vary between 55~225 mm Hg, to choose required value use buttons  or  .

Current Date and Time Setting


Fig. 30. Time settings. Press corresponding button to set system time.

Start/Finish of therapy cycle

Press  to start therapy. Press  to finish treatment. Sensor display will be automatically blocked after 1 min. therapy interval.


Fig. 31. Therapy in continuous mode.


Fig. 32. System state: therapy in continuous mode.

4 h: time passed since therapy started, display shows 0 h if therapy time is set as endless (∞ symbol).

8 h: remained therapy time, "Not applicable" (N/A) means that therapy time is set on "non limit" option, the device will work until user turn off the device, until battery charge is over, until automatic turn off of the device.

125 mm Hg: the set pressure value is displayed in continuous therapy mode. The value may vary between 55~225 mm Hg.

0 mm Hg: common pressure value is displayed in continuous therapy mode.


Fig.33. System state: Intermittent therapy mode.

5 min: therapy interval at high pressure in intermittent mode.

The value can vary within 2 min\3 min\4 min\ 5 min.

2 min: therapy interval at low pressure in intermittent mode. Available values to choose are 1 min\2 min\ 3 min.

125 mm Hg: the set pressure in intermittent mode. Value can be set between 55 mm Hg ~ 225 mm Hg.

0 mm рт.ст.: common pressure value in intermittent therapy mode is displayed.

Current state "Unblocked". Press the button to block sensor display as the therapy is on.

Current state "Blocked". Press the button to unblock sensor display, confirming message will appear and corresponding image.


Fig. 34. Unblock of operation screen.

Press , then to confirm sensor display unblock.

Device Disconnection

Attention! The device can be disconnected from drainage tubing for a short period for bathing.

Attention! According to clinical instructions change dressing kit if treatment is interrupted or stopped for 2 hours or more.

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1. To stop therapy cycle press button on sensor display.
2. Press button to switch off the device, disconnect it from the main power.
3. Продвиньте зажим трубки в направлении коннектора трубки.
4. Move the tube clamp in direction of tube connector.
5. Close the tube clamp tightly to prevent contents leakage into drainage tube: several clicks must be heard.
6. Pull the connector aside to disconnect drainage tube from the canister.
7. Close the drainage tube end with gauze material to collect liquid, flowing out of drainage tube.

Information Review

Information System Interface

Contents	Operation instructions
Recent therapy records	Press  → RECORD for data review.
Device components operation data	Press  → DEVICE for data review..
Contact data Version	Press  → VERSION for data review.

Exudate Information

Press the chart "Exudate volume" to move to corresponding window.


Fig. 36. Window of exudate volume data.

-  Pressing this button allows to check information on exudate volume.
-  Pressing this button allows to delete current records on exudate volume, before records deletion there is message asking to confirm the action.

9. Troubleshooting

Warnings and Alarm Signals

Warnings and alarms	Cause	Remedy
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Low Bat" appears in the middle of display.	Low battery	Urgent! Charge the device connecting it to the mains power.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "System Error" appears in the middle of display.	1. System error. 2. Current vacuum pressure value is much higher than the set vacuum pressure value.	1. Check gas circuit connection. 2. Wait until the vacuum pressure drops to zero, then restart the device. 3. Check and ensure that tubing is not kinked, crimped or blocked in any way. Wait until the vacuum pressure drops to zero and restart the device.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Canister Full" appears in the middle of display.	Canister is overfilled	Urgent! Replace the canister for a new one, no later than in 30 min*.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Time's Up" appears in the middle of display.	Treatment time is up.	Press MENU for a while to disable the alarm and skip to the main menu.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Too Hot" appears in the middle of display.	The system temperature is too high.	Device must be turned off for at least 30 minutes before restart to prevent permanent damage to device.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Too Cold" appears in the middle of display.	The system temperature is too low.	Place the device in a warmer environment.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "Sealing Error" appears in the middle of display.	Sealing of wound dressing or air connections is seriously (high vacuum leakage).	Check and ensure that the wound bed, canister and tubing are well sealed.
Red LED flashes on the front panel and audible alarm sounds each 0,5 sec. The message "Alarm!" appears in the display right bottom angle. The message "No Operation" appears in the middle of display.	The device will enter standby mode when no operation for too long time.	When alarm occurs, press any button to disable the alarm.
Yellow LED flashes on the front panel each 0,5 sec. and audible alarm sounds each 20 sec. The message "Charge" appears in the display right angle.	Low battery.	Plug device into AC power to start charging immediately.
Yellow LED flashes on the front panel each 0,5 sec. and audible alarm sounds each 20 sec. The message "Sealing Error" appears in the display right angle.	Sealing of wound dressing or air connections is seriously (minor vacuum leakage).	Check and ensure that the wound bed, canister and tubing are well sealed.
Yellow LED constantly lights on the front panel. The message "Power Error" appears in the bottom right angle. Symbol appears in the bottom left angle of the display.	The AC adapter input voltage is abnormal.	Unplug the AC adapter and check the voltage of AC power supply. The AC adapter is broken. Replace it with a new one that has the same specifications.
The message "No Battery" appears in the bottom right angle of the display.	No battery or battery invalidation.	Check whether the battery is connected. Replace the battery.
No feedback on button pressing. Symbol appears in the bottom left angle of the display.	Keypad is blocked.	To unblock the keypad press Menu, and then Lock.

*Use canister intended especially for Sensivac device.

Troubleshooting

Symptom	Possible cause	Remedy
Unable to start the device	Battery Empty	Recharge immediately.
Auto power off after startup	Low battery	Recharge immediately.
	Long key press of the power button	Once the device enters into the startup screen, release the power button and ensure that the button is not stuck.
Unable to monitor exudate volume	Inappropriate canister installation	Ensure that there is no clearance between canister and the device; Push the canister firmly into the device.
	The device is not flat after placement.	Place the unit on a flat horizontal surface.

Attention! Do not try to troubleshoot at your own discretion if troubleshooting measures above do not work or other failures occur. Please contact your supplier immediately.

10. Service

Attention! Do not attempt to open the enclosure. Please Contact your representative, distributor or an authorized provider if service is required.

Dressing kits and canisters used with the device are intended for single use only and do not require any technical service. All disposable items (used canister kits and used dressings) should be handled in accordance with local medical waste disposal regulations.

- Wipe the device surface with soft wet cloth wetted with non-abrasive washing solution, avoiding direct contact with liquids especially close to the holes on the right side of the device and vacuum port. Wipe the device with dry cloth;
- Assure that the device is stored clean and dry.

Attention! No condensed moisture or liquid should drip into the device. Do not immerse the device in water or expose it to water sprays or other liquids.

When the device is returned to service centre for repair it must be preliminary cleaned and disinfected. The device and its components are considered dangerous and potentially infected. The device must be properly marked with special label and be stored in isolated room before disinfection.

Use personnel protection equipment while device cleaning to prevent risk of infection and reduce contact with blood and other biological liquids.

Note: The device contains accumulator battery which is a power source for system settings (time, date) when the device is off. If the batteries are replaced, corresponding time and date should be set.

Battery Maintenance

The device contains a lithium polymer rechargeable battery. The device can run on battery power and/or while plugged in with the power adapter. During therapy, if the battery indicator shows that the battery is low, plug device into AC power to start charging. The device is designed to permit use of the device while the internal battery is charging. Charging will automatically stop once the battery is fully charged.

If the device is fully charged and is not going to be used, unplug the AC adapter from line power.

Attention! For the first time the battery should be charged for over 4 hours. Charge the battery completely each time before use.

11. Transportation and storage

The device contains precision electronic components, so avoid vibration during transport. Keep the device at a normal temperature and avoid prolonged sun exposure during transport.

During transportation and storage, following climatic conditions should be observed:

- Storage temperature: from -10 °C to + 55°C
- Relative air humidity up to 95 %
- Atmosphere pressure 500 gPa to 1060 gPa

When the suction device is stored and transported under low temperature conditions, leave it before operation for 12 hours under room temperature to avoid the damage of the suction device pump.

Transportation of aspirators for vacuum wound therapy SENSITEC should be performed by roofed transport of any type according to shipping rules relative to using type of transport. Transportation by rail should be made by container shipments.

12. Disposal

The suction device must be disposed separately from household waste for environmental protective purpose. The suction device must be disposed in compliance with Directive EC 2002/96/EC – WEEE (Waste Electrical and Electronic Equipment). Please, address all questions, if any, to Local Utility Service responsible for the waste disposal.

13. Accessories/Set Contents

1. Canister (collection bowl) with gelling agent and cups for holes packing, volume 350 ml.
2. Canister (collection bowl) with gelling agent, volume 800 ml.
3. Holding tool for canister of the volume 800 ml.
4. Canister's cover (for volume 800 ml) with cups for holes packing.
5. Bacteria filter.
6. Connecting pipe with end cups.
7. Connector, angular.
8. Bracket for rail fastening.
9. Wall rail.
10. Transfer bag.

Unpack the device accurately and assure that all enumerated components are included. If any component is absent, appeal to Supplier for support.

14. EMC Declaration

Guidance and manufacturer's declaration-electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that it is used in such environment.		
Emission test	Standard	Electromagnetic environment – guidance
RF emissions CISPR 11	Group 1	The pump only uses radio frequency energy when running internal functions. The pump's radio frequency emissions are therefore very low and will not cause any electromagnetic interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	
Harmonic emissions IEC61000-3-2	Class A	
Voltage fluctuations and scintillation IEC 61000-3-3	Compliant	

Guidance and statements regarding electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such environment.

Immunity test	IEC60601 test level	Standard	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 Kv contact discharge ±15 kV air discharge	±8 kV contact discharge ±15 kV air discharge	Flooring must be wood, concrete or ceramic tile. If the floor is lined with synthetic materials there must be a relative humidity of at least 30%.
Electrical fast transient (EFT) IEC 61000-4-4	±2 kV power cord ±1 kV I/O cable	±2 kV power cord	The network power source must be of typical commercial or hospital quality.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	
Voltage drops, short interruptions and changes IEC 61000-4-11	<5% UT(drop >95% UT) 0.5 cycles 40% U (drop 60% UT) 5 cycles 70% UT(drop 30% UT) 25 cycles <5% UT (drop >95% UT) 5 seconds	<5% UT(drop >95% UT) 0.5 cycles 40% U (drop 60% UT) 5 cycles 70% UT(drop 30% UT) 25 cycles <5% UT(drop >95% UT) 5 seconds	The network power source must be of typical commercial or hospital quality. If the pump needs to run continuously we recommend using an uninterruptible power supply (UPS) in case of interruptions to mains power.
Power frequency magnetic field (50/60 Hz) IEC 61000-4-8	400 A/m	400 A/m	The power frequency magnetic field must be at a typical level for typical commercial or hospital environments.

Note: U_T refers to the voltage of the AC power network before voltage testing

Guidance and statements regarding electromagnetic immunity

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such environment.

Immunity test	IEC 60601 Test level	Standard	Electromagnetic environment – guidance
Conducted immunity IEC61000-4-6	3 Vrms 150 k–80 MHz	3 V	Portable and mobile radio frequency communications devices must be used at the stipulated distance away from the equipment or any system components (including the cables). This isolation distance is calculated using appropriate formulas selected on the basis of the radiation frequency. Recommended calculation formulas for isolation distance are: $d=1.2\sqrt{P}$ $d=1.2\sqrt{P}$ 80 M–800 MHz $d=2.3\sqrt{P}$ 800 M–2.5 GHz where P is the rated maximum output power of the transmitter, in W. d is the recommended distance, in m. The field strength of the radio frequency transmitter obtained by measuring the electromagnetic field must be within compliant levels for every frequency range ^a . May cause interference if used in close proximity to equipment with the following symbols: 
Radiated immunity IEC61000-4-3	3 V/m 80 M–2.5 GHz	3 V/m	

Note 1: Between 80 MHz–800 MHz, use a formula for higher frequency bands.

Note 2: The above guidance is not suitable for use in all conditions. Material structures, objects and people can absorb and reflect electromagnetic waves, affecting electromagnetic propagation.

^a Field strengths for radio (honeycomb and wireless) handset base stations and terrestrial mobile radio receiving apparatus, antenna reception apparatus, and FM and AM radio/television broadcasts cannot be accurately estimated using a purely theoretical approach.

The use of electromagnetic field measuring methods should be considered when estimating the electromagnetic environment produced by a fixed radio frequency transmitter. If the measured field strengths in the environment in which this pump is used exceed stipulated RF levels, the pump must be observed to check whether normal operation is possible. Relevant action must be taken as soon as abnormalities are discovered, such as changing the position of the pump or moving it to another environment.

^b In a frequency range of between 150 k–80 MHz, field strength should be less than 3 V/m.

It is recommended that the pump be kept away from portable/mobile RF communications equipment.

The device is intended for use in the electromagnetic environment specified below. The customer or the user of the device should assure that they are used in such environment.

Transmitter ratings Maximum output power (W)	Calculate isolation distances on the basis of transmitter frequencies (m)		
	150 k–80 MHz $d=1.16\sqrt{P}$	80 M–800 MHz $d=1.16\sqrt{P}$	800 M–2.5 GHz $d=2.33\sqrt{P}$
0,01	0,12	0,12	0,23
0,1	0,38	0,38	0,73
1	1,2	1,2	2,3
10	3,8	3,8	7,3
100	12	12	23

If the rated maximum output power of the transmitter is not included above, the isolation distance may be estimated by using the formula in the corresponding column. In the formula, P is the rated maximum output power of the transmitter as provided by the manufacturer, in W.

Note 1: Between 80M-800MHz, use a formula for higher frequency bands.

Note 2: The above guidance is not suitable for use in all conditions. Material structures, objects and people can absorb and reflect electromagnetic waves, affecting electromagnetic propagation.

15. Limited Warranty

Warranty term for the device is 12 month from purchase date, proved by stamp and signature of the dealer in the device passport. The warranty must be accompanied by the specified documents.

Contact Information

Apexmed International BV (Keizersgracht 62-64, 1015 CS, Amsterdam, the Netherlands | w: apexmed.nl | e: info@apexmed.nl | t: +31 20 520 74 03 | f: +31 20 520 75 10).

■ Foryou Medical Electronics Co., Ltd.(No.1, North Shangxia Road, Dongjiang Hi-tech Dist., Huizhou City, Guangdong, 516000, China | w: ForyouNPWT.com | e: info@ForyouNPWT.com | t: +86-752-5300546 | f: +86-752-5300666)

The company, authorized by manufacturer to accept and satisfy customer's requirements regarding non-conforming products: "MEDICOM Co., Ltd.", Russian Federation 143570, Moscow region, Istra district, Novopetrovskoe settlement, pervomayskaya str., 59, Tel./Fax: (495) 651-9116, e-mail: apexmed@apexmed.ru

All rights reserved.

Partial or complete reproduction without written permission of the rightholder is forbidden.

© 2004 Apexmed

© 2006 Sensitec

16. Symbols

Warning: Misuse could cause personal injury or property damage

Useful information

Caution! Appeal to accompanying documentation

B-type equipment

Temperature limits

Caution! Laser emission

Class II device

Disposal according to the set rules

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