CALIFORNIA ALL-PURPOSE ACKNOWLEDGMENT

CIVIL CODE § 1189

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California)
County of Ventura)

On March 7, 2023 before me, Shacoya Braggs, Notary Public
Date Here Insert Name and Title of the Officer
personally appeared Dolphine Amanya Katiwa, PhD
Name(s) of Signer(s)

who proved to me on the basis of satisfactory evidence to be the person(s) whose name(s) is/are subscribed to the within instrument and acknowledged to me that he/she/they executed the same in his/her/their authorized capacity(ies), and that by his/her/their signature(s) on the instrument the person(s), or the entity upon behalf of which the person(s) acted, executed the instrument.

I certify under PENALTY OF PERJURY under the laws of the State of California that the foregoing paragraph is true and correct.

WITNESS my hand and official seal.

Signature Shacoya Braggs
Signature of Notary Public

Place Notary Seal Above

OPTIONAL

Though this section is optional, completing this information can deter alteration of the document or fraudulent reattachment of this form to an unintended document.

Description of Attached Document

Title or Type of Document: Attestation Statement Document Date: March 7, 2023
Number of Pages: _____ Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____
☐ Corporate Officer — Title(s): _____
☐ Partner — ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other: _____
Signer Is Representing: _____

Signer's Name: _____
☐ Corporate Officer — Title(s): _____
☐ Partner — ☐ Limited ☐ General
☐ Individual ☐ Attorney in Fact
☐ Trustee ☐ Guardian or Conservator
☐ Other: _____
Signer Is Representing: _____

Medtronic

18000 Devonshire Street
Northridge, CA 91325 USA
www.medtronic.com
Tel 818-576-5211

ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the following Instruction For Use, attached herein, to be a true and accurate copies.

Dolphy Katiwa

Dolphine Amenya Katiwa, PhD
Senior Regulatory Affairs Manager

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
 Bluetooth®

Medtronic

Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("TM") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.

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Medtronic Australasia Pty. Ltd.
Tel: 1800 668 670

Bangladesh:

Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-9903995417
or (+880)-1714217131

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
Guardian 4

Introduction

The Guardian 4 transmitter (MMT-7841Q) with Bluetooth® wireless technology is a component of the continuous glucose monitoring (CGM) system with Bluetooth® wireless technology.

Indications for use

The Guardian 4 transmitter (MMT-7841Q) is a rechargeable device that powers the glucose sensor. The transmitter collects and calculates sensor data and sends the data via Bluetooth® wireless technology to a compatible display device, such as a mobile application, for the management of diabetes mellitus. The transmitter is only compatible with the Guardian 4 sensor (MMT-7040Q) and is indicated for single-patient or multiple-patient use. The transmitter is intended for persons ages 7 years and older.

Contraindications

No contraindications are associated with Guardian 4 transmitter use. For contraindications related to CGM, see the Guardian App User Guide.

Clinical Benefits

The Guardian 4 transmitter is a component of the CGM system that provides sensor glucose values. See the Guardian App User Guide for the clinical benefits of systems that use the Guardian 4 transmitter.

User Safety

Warnings

- Always refer to the Guardian 4 Sensor User Guide for all precautions, warnings, and instructions related to the sensor. Not referring to the Guardian 4 Sensor User Guide can result in serious injury or damage to the sensor.
- Do not allow children to put small parts in their mouth. This product may pose a choking hazard that can result in serious injury or death.
- Do not change or modify the device unless expressly approved by Medtronic Diabetes. Modifying the device can cause serious injury, interfere with the ability to operate the device, and void the warranty.
- Do not expose the transmitter to Magnetic Resonance Imaging (MRI) equipment, diathermy devices, or other devices that generate strong magnetic fields (for example x-ray, CT scan or other types of radiation). Exposure to a strong magnetic field has not been evaluated and can cause the device to malfunction, result in serious injury, or be unsafe. If the transmitter is exposed to a strong

magnetic field, discontinue use and contact a local Medtronic support representative for further assistance.

- Do not use the tester if it comes in contact with blood. Touching blood can cause infection.
- Bleeding may occur after inserting the sensor. Always make sure that the site is not bleeding before connecting the transmitter to the sensor. Blood can get into the transmitter connector and damage the device. Discard the device if damaged. If bleeding occurs, apply steady pressure with a sterile gauze, pad, or clean cloth at the insertion site until bleeding stops. After bleeding stops, connect the transmitter to the sensor.
- Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite and result in serious injury.
- If a serious incident related to the device occurs, immediately report the incident to the manufacturer and local competent authority.
- For questions or concerns related to product use, contact a local Medtronic support representative for assistance.
- For medical questions or concerns, contact a healthcare provider.

Precautions


- Do not use the transmitter adjacent to other electrical equipment that may cause interference with the normal system operation.
- Only use the Guardian 4 sensor (MMT-7040Q) with the transmitter. Do not use any other sensor. Other sensors are not intended for use with the transmitter and will damage the transmitter and the sensor.
- Only use the green colored tester (MMT-7736L) with the transmitter. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.

Figure 1. Transmitter and tester

- Always use the tester when cleaning the transmitter. Do not use any other test plug with the transmitter. Use of another test plug can allow water to get into the transmitter or can prevent proper cleaning. Water can damage the transmitter.
- Do not twist the tester or sensor while attached to the transmitter. Twisting the tester or sensor will damage the transmitter.
- Do not allow the tester to come in contact with any liquid when not connected to the transmitter. A wet tester can damage the transmitter.
- Do not allow the transmitter to come in contact with any liquid when not connected to a sensor or to the tester. Moisture will damage the transmitter and a wet transmitter can damage the sensor.
- Do not clean the O-rings on the tester with any substances. Cleaning the O-rings can damage the tester.

Figure 2. O-rings

IEC 60601-1-2:2014, 4th Edition; Special EMC Precautions for Medical Electrical Equipment

1. Special Precautions regarding Electromagnetic Compatibility (EMC): This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment where common levels of radiated "E" (V/m) or "H" fields (A/m) exist, such as cellular phones, Wi-Fi™, Bluetooth® wireless technology, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
2. Portable and mobile RF communications equipment can affect medical electrical equipment. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is causing the interference.

3. Be careful when using the transmitter closer than 30 cm (12 in) to portable radio frequency (RF) equipment or electrical equipment. If the transmitter must be used next to portable RF equipment or electrical equipment, observe the transmitter to verify correct system operation. Degradation of the performance of the transmitter could result.
4. The essential performance (EP) of the transmitter is to measure and transmit to a monitoring device the sensing device's signal value(s) within the transmitter's accuracy requirements under the specified use conditions outlined in the Guardian App User Guide and for the duration of the expected service life. If the transmitter experiences electromagnetic disturbances, either no or incorrect data may be transmitted. In such situations, refer to the operation, maintenance, and troubleshooting instructions within the applicable user guides. You may also use the tester to test if the transmitter is operating properly. If the transmitter is damaged or if it cannot communicate with the display device, contact your local Medtronic support representative for assistance.

Assistance

Contact a local Medtronic support representative for assistance and if a copy of the Guardian App User Guide is needed.

Risks and side effects

The product contains small parts and may pose a choking hazard that can result in serious injury or death.

Side effects may include discomfort and skin irritation at the insertion site.

Hazardous substances

None.

Allergens

None known.

Using the transmitter

Components needed

- Guardian 4 transmitter (MMT-7841Q)
- Tester (MMT-7736L)
- Charger (MMT-7715)

Figure 3. Components

Preparing the transmitter


The transmitter contains a non-replaceable, rechargeable battery that can recharge as needed with the charger. The transmitter needs to be charged before use. The charger has a green light that shows the charging status and a red light that communicates any problems during charging. If there is a red light, see *Troubleshooting*, page 31. The charger requires one AAA alkaline battery.

Note: If the battery is installed incorrectly or is low, the charger will not work. Repeat the battery installation steps using a new battery.

Installing a battery in the charger

To install a battery in the charger:

1. Push the battery cover in and slide it off (as shown in image A in step 3).
2. Insert a new AAA alkaline battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger.
3. Slide the cover back on the charger until it clicks into place (as shown in image B in step 3).


Charging the transmitter

CAUTION: Always charge the transmitter before inserting the sensor. A depleted transmitter does not function. A fully charged transmitter works at least seven days without recharging. A depleted transmitter can take up to two hours to recharge.


CAUTION: Do not store the transmitter on the charger for more than 60 days. Disconnect and reconnect to the charger to re-charge again before use. If the transmitter is left on the charger for more than 60 days, the transmitter battery will be permanently damaged.

To charge the transmitter:

1. Push the transmitter and the charger together to connect the transmitter to the charger.

2. Within 10 seconds after the transmitter is connected, a green light on the charger will flash for one to two seconds as the charger powers on. For the rest of the charging time, the green light on the charger will continue to flash in a pattern of four flashes with a pause between the four flashes.
3. When charging is complete, the green light on the charger will stay on, without flashing, for 15 to 20 seconds and then turn off.

4. After the green charger light turns off, disconnect the transmitter from the charger. The green light on the transmitter starts to flash.

Pairing the transmitter

The transmitter must be paired to the system before a sensor can be used. Always refer to the Guardian App User Guide for instructions on how to pair the transmitter to the system.

Inserting the sensor


Always refer to the Guardian 4 Sensor User Guide for instructions on how to insert the sensor.

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
Before proceeding, have the Guardian App User Guide available.

To connect the transmitter to the sensor:

1. After the sensor is inserted, consult the Guardian 4 Sensor User Guide for details on how to apply the required tape before connecting the transmitter.
2. Hold the rounded end of the inserted sensor to prevent it from moving during connection.

3. Hold the transmitter as shown. Line up the two notches on the transmitter with the side arms of the sensor. The flat side of the transmitter should face the skin.

4. Slide the transmitter onto the sensor connector until the sensor arms snap into the notches on the transmitter. If the transmitter is properly connected, and if the sensor has had enough time to become hydrated with the interstitial fluid, the green light on the transmitter will flash 6 times.


Note: If the transmitter does not flash, see *Troubleshooting*, page 31.

5. When the transmitter light flashes green after connecting to the sensor, use the system to start the sensor. For more instructions, see the Guardian App User Guide.
6. Attach the adhesive tab of the sensor to the transmitter.
7. After the transmitter is connected, consult the Guardian 4 Sensor User Guide for instructions on how to apply a second tape.
8. Refer to the Guardian App User Guide for pairing.

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To disconnect the transmitter from the sensor:

1. Carefully remove any tape from the transmitter and sensor.

2. Remove the adhesive tab from the top of the transmitter.
3. Hold the transmitter as shown, and pinch the flexible side arms of the sensor between the thumb and forefinger.
4. Gently pull the transmitter away from the sensor.

Removing the sensor

Always refer to the Guardian 4 Sensor User Guide for instructions on how to remove the sensor.

Reconnecting the transmitter to a sensor that is already inserted

The transmitter can be reconnected to the sensor currently in use. Simply connect the transmitter to the sensor that is already inserted. Confirm Reconnect Sensor when the display device detects the transmitter. It may take a few seconds to establish a connection when reconnecting a sensor. Reattach the adhesive tab of the sensor to the transmitter and reapply any required tape. When reconnected, the sensor goes through another warm-up period.


Tester

The tester is used to test the transmitter to make sure it is working. The tester is also used as a required component to create a waterproof seal when cleaning the transmitter. Properly connecting the tester to the transmitter ensures that fluids do not come in contact with the connector pins inside the transmitter. Fluids can cause connector pins to corrode and affect the performance of the transmitter.

Do not twist the tester while attached to the transmitter. This will damage the transmitter.


The tester can be used for one year. If the tester is used for more than one year, the connector pins inside the transmitter can be damaged, because the tester cannot continue to provide a waterproof seal. For instructions on how to check the connector pins, see *Inspecting the transmitter connector pins*, page 20.

CAUTION: Only use the green colored tester (MMT-7736L) with the transmitter. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.

Inspecting the transmitter connector pins

This image is an example of how the connector pins should look for the transmitter.

Look inside the connector opening of the transmitter to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger or the display device. Contact a local Medtronic support representative. It may be time to replace your transmitter.


Look for moisture inside the connector opening. If any moisture is present, allow the transmitter to dry for at least one hour. Moisture inside the connector opening could cause the transmitter to not work properly and could cause corrosion and damage over time.

Connecting the tester for testing or cleaning

Before proceeding, have the Guardian App User Guide available.

To connect the tester:

1. Hold the transmitter body as shown and pinch the side arms of the tester.


2. Push the tester into the transmitter until the flexible side arms of the tester click into the notches on both sides of the transmitter.

When properly connected, the green light on the transmitter flashes 6 times.

3. To test the transmitter, check the sensor icon in the app to ensure that the transmitter is sending a signal (see the Guardian App User Guide).
4. To clean the transmitter, see *Cleaning the transmitter*, page 22.
5. After testing or cleaning, disconnect the tester from the transmitter.

Disconnecting the tester**To disconnect the tester:**

1. Hold the transmitter body as shown and pinch the side arms of the tester.

2. With the tester arms pinched, gently pull the transmitter away from the tester.

Note: To save transmitter battery life, do NOT leave the tester connected after cleaning or testing.

Cleaning the transmitter

The transmitter is intended for personal use at home (single-patient use) or for use in healthcare facilities (multiple-patient use). Single-patient use requires cleaning after each use, while multiple-patient use requires cleaning and disinfection after each use. When using the transmitter in a healthcare facility, always follow the cleaning and disinfecting procedure for multiple-patient use.

WARNING: Do not discard the transmitter in a medical waste container or expose it to extreme heat. The transmitter contains a battery that may ignite and result in serious injury.

Note: The tester is a required component for cleaning the transmitter. For details, see *Tester*, page 19.

CAUTION: Do not use an automated washer-disinfector to clean or disinfect the device. Using an automated washer-disinfector to clean or disinfect the device will cause damage to the transmitter.

For single-patient use

Always clean the transmitter after each use.

To clean the transmitter, use these materials:

- mild liquid soap
- soft-bristled toddler toothbrush
- container
- clean, lint-free dry cloths


Use life

The transmitter can be cleaned up to 122 times or for one year, whichever comes first. Discard the transmitter at this point. If the transmitter is used beyond 122 times or one year, the cleaning process may damage the device. Contact a local Medtronic support representative to order a new transmitter.


WARNING: Do not use the device if there is any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter and result in serious injury. Call a local Medtronic support representative and discard the device according to local regulations for battery disposal (non-incineration).

To clean the transmitter:

1. Wash hands thoroughly.
2. Attach the tester to the transmitter to create a waterproof seal.


3. If there is adhesive residue on the transmitter, see *Removing adhesive residue*, page 30.
4. Rinse the transmitter under room temperature tap water for at least one minute, and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.


5. Prepare a mild liquid soap solution using 5 mL (1 teaspoon) of mild liquid soap per 3.8 L (1 gallon) of room temperature tap water.


6. With the tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.


7. Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.

8. Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.


9. Dry the transmitter and tester with a clean, dry cloth.

10. Place the transmitter and tester on a clean, dry cloth and air dry them completely.

11. Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

For multiple-patient use

When using the transmitter in a healthcare facility, always clean and disinfect the transmitter after each use.

WARNING: Standard precautions must be followed when handling or using this device. All parts of the system should be considered potentially infectious and are capable of transmitting blood-borne pathogens between patients and healthcare professionals.

The transmitter must be disinfected after use on each patient. This system may only be used for testing multiple patients when Standard Precautions and disinfection procedures provided by Medtronic Diabetes are followed.

To clean the transmitter, use these materials:

- gloves
- mild liquid soap

- soft-bristled toddler toothbrush
- 8.25% bleach
- two containers
- clean and lint-free dry cloths

Use life


The transmitter can be cleaned and disinfected up to 122 times or for one year, whichever comes first. Discard the transmitter at this point. If the transmitter is used beyond 122 times or one year, the cleaning and disinfection process may damage the device. Contact a local Medtronic support representative to order a new transmitter.

To clean and disinfect the transmitter:


1. Wash hands and put on gloves.
2. Inspect the inside of the transmitter's connector opening for any sign of body fluid. For instructions on how to inspect the connector pins, see *Inspecting the transmitter connector pins*, page 20.

CAUTION: The person inspecting the transmitter must have sufficient vision that enables him or her to see small drops of body fluid or debris.

WARNING: If there is any body fluid in the connector opening, the transmitter must be discarded. Because the transmitter contains a battery, do not discard in a bio-waste container. Instead, continue to clean and disinfect the transmitter, and then discard according to local regulations for battery disposal (non-incineration).


3. Attach the tester to the transmitter to create a waterproof seal.

4. If there is adhesive residue on the transmitter, see *Removing adhesive residue*, page 30.
5. Rinse the transmitter under room temperature tap water for at least one minute and until visibly clean. Make sure all hard-to-reach areas are rinsed completely.


6. Prepare a mild liquid soap solution using 5 mL (1 teaspoon) of mild liquid soap per 3.8 L (1 gallon) of room temperature tap water. Make sure to prepare a fresh solution for each use.
7. With the tester still attached, submerge the transmitter in the mild liquid soap solution and soak for one minute.


8. Holding the tester, brush the entire surface of the transmitter using a soft-bristled toddler toothbrush. Make sure to brush all hard-to-reach areas until visibly clean.


9. Rinse the transmitter under running room temperature tap water for at least one minute, and until all visible liquid soap is gone.


10. Dry the transmitter and tester with a clean, dry cloth.

11. Prepare a 1:10 bleach solution by using one (1) part 8.25% bleach to nine (9) parts water, for a final concentration of 0.8%. Make sure to prepare a fresh solution for each use.
12. Make sure that the previous cleaning steps are completed before disinfection. With the tester still attached, soak the transmitter in the bleach solution for 20 minutes.


13. Rinse the transmitter under running room temperature tap water for three minutes.

14. Place the transmitter and tester on a clean, dry cloth and air dry them completely.

WARNING: If there is any body fluid inside the connector opening on earlier inspection, the transmitter must be discarded with the tester still attached, according to local regulations for battery disposal (non-incineration).

15. Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

16. Inspect the housing of the transmitter for any signs of cracking, flaking, or damage. If these signs are present, the transmitter must be discarded according to local regulations for battery disposal (non-incineration).

WARNING: Do not use the device if there is any cracking, flaking, or damage to the housing. Cracking, flaking, or damage to the housing are signs of deterioration. Deterioration of the housing can affect the ability to properly clean the transmitter, and result in serious injury. Call a local Medtronic support representative and discard the device according to local regulations for battery disposal (non-incineration).

17. Discard the used gloves and thoroughly wash hands with soap and water.

Removing adhesive residue


Follow these instructions if there is adhesive residue present on the transmitter.

Use cotton swabs and a medical adhesive remover such as Detachol™*, a mineral spirit, to remove adhesive residue.

Note: During testing, Detachol™* was used to remove the adhesive residue from the transmitter. Detachol™* is recommended for use but may not be available in all countries.

To remove adhesive residue:

1. Make sure the tester is attached to the transmitter.
2. Soak a cotton swab in the medical adhesive remover.
3. Hold the tester and gently rub the adhesive remover on the transmitter until the residue is removed.

4. Continue with the cleaning procedure. See *Cleaning the transmitter*, page 22 for details.

Cleaning the charger

This procedure is for general cleaning as required, based on physical appearance.

CAUTION: Do not immerse the charger in water or any other cleaning agent. The charger is not waterproof. Water can damage the charger and cause the device to malfunction.

To clean the charger:

1. Wash hands thoroughly.
2. Use a damp cloth with mild cleaning solution, such as a dishwashing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents, such as paint thinner or acetone, to clean the charger.
3. Place the charger on a clean, dry cloth and air dry for two to three minutes.

Bathing and swimming

After the transmitter and sensor are connected, they form a waterproof seal to a depth of 2.4 meters (8 feet) for up to 30 minutes. Shower and swim without removing them.

Troubleshooting

The table shown contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see the Guardian App User Guide.

Table 1. Troubleshooting issues

Problem	Likely Cause(s)	Resolution
The transmitter is connected to the charger and no lights come on.	The transmitter connector pins are damaged or corroded. The charger battery has no power or no battery is inserted.	1. Check the transmitter connector pins for damage or corrosion. For more information about the connector pins, see <i>Inspecting the transmitter connector pins, page 20</i>. If the pins are damaged or corroded, contact a local Medtronic support representative. It may be time to replace the transmitter. 2. If there is no damage to the connector pins, replace the battery in the charger. For instructions on replacing the charger battery, see <i>Installing a battery in the charger, page 15</i>.
During charging, the flashing green light on the charger turns off and a longer flashing red light appears on the charger.	The charger battery is low on power.	Replace the battery in the charger. For instructions on replacing the charger battery, see <i>Installing a battery in the charger, page 15</i> .
When charging, the flashing green light on the charger turns off and there is a series of quick flashing red lights on the charger for two seconds at a time.	The transmitter is low on power.	1. Charge the transmitter continuously for one hour. If flashing does not stop, proceed to step 2. 2. Charge the transmitter continuously for eight hours. If flashing does not stop, call a local Medtronic support representative. It may be time to replace the transmitter.
During charging, a mix of quick and long flashing red lights appear on the charger.	The charger and the transmitter are low on power.	1. Replace the battery in the charger. For instructions on replacing the charger battery, see <i>Installing a battery in the charger, page 15</i>. 2. Charge the transmitter continuously for one hour. If the quick flashing red lights do not stop, proceed to step 3. 3. Charge the transmitter continuously for eight hours. If flashing does not stop, call a local Medtronic support representative.

Table 1. Troubleshooting issues (continued)

Problem	Likely Cause(s)	Resolution
		tative. It may be time to replace the transmitter.
When connected to the sensor, the green light on the transmitter does not flash.	The transmitter is not fully connected. The transmitter is low on power. The sensor is not properly inserted into the body.	1. Disconnect the transmitter from the sensor. 2. Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3. Fully charge the transmitter and connect it to the tester. If the green light still does not flash, see troubleshooting on "When connected to the tester, the green light on the transmitter does not flash." If the green light flashes, proceed to step 4. 4. Disconnect the transmitter from the tester, wait at least five seconds, and connect the transmitter to the sensor. If the green light still does not flash, proceed to step 5. 5. The sensor may not be properly inserted into the body. Remove the sensor from the body and insert a new sensor.
When connected to the tester, the green light on the transmitter does not flash.	The transmitter is not fully connected. The transmitter is low on power.	1. Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2. Fully charge the transmitter. 3. Test the transmitter with the tester again. If there is no flashing green light, call a local Medtronic support representative. It may be time to replace the transmitter.
The transmitter battery does not last for seven days.	The transmitter is not fully charged when connected to the sensor. The transmitter and display device frequent-	1. Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for the duration of one sensor use, proceed to step 2. 2. Move away from any equipment that can cause RF interference. For more informa-

Table 1. Troubleshooting issues (continued)

Problem	Likely Cause(s)	Resolution
	ly lose wireless connection.	tion on RF interference, see the Radio Compliance Information sheet included with the display device. 3. Make sure the display device and the transmitter are located on the same side of the body to minimize any RF interference. If a fully charged transmitter battery continues to lose power before a full seven days, call a local Medtronic support representative. It may be time to replace the transmitter.
The transmitter loses connection with the display device.	The display device is out of range. There is RF interference from other devices.	1. Move away from any equipment that can cause RF interference. For more information on RF interference, see the Radio Compliance information sheet included with the transmitter. If the transmitter is still not communicating with the display device, proceed to step 2. 2. Make sure the display device and the transmitter are located on the same side of the body to minimize any RF interference. If the transmitter is still not communicating with the display device, call a local Medtronic support representative for assistance.
Note: An alarm or alert occurs and a message appears when the transmitter has lost connection with the display device for 30 minutes.		

Storage

Store the transmitter, charger, and tester in a clean, dry location at room temperature. If the transmitter is not in use, you must charge the transmitter at least once every 60 days.

CAUTION: Do not store the transmitter on the charger. If the transmitter is left on the charger for more than 60 days, the battery will be permanently damaged.

Disposal

Do not dispose of the transmitter, charger, and tester in unsorted municipal waste stream. Dispose of the transmitter, charger, and tester according to local regulations for electronic waste disposal.

Technical Specifications

Table 2. Product specifications

Biocompatibility	Transmitter: Complies with EN ISO 10993-1
Applied parts	Transmitter Sensor
Operating conditions	Transmitter temperature: 0°C to 45°C (32°F to 113°F) Caution: When operating the transmitter on a tester in air temperatures greater than 41°C (106°F), the temperature of the transmitter may exceed 43°C (109°F). Transmitter relative humidity: 10% to 95% with no condensation Transmitter pressure: 57.60 kPa to 106.17 kPa (8.4 psi to 15.4 psi) Charger temperature: 10°C to 40°C (50°F to 104°F) Charger relative humidity: 30% to 75% with no condensation
Storage conditions	Transmitter temperature: -20°C to 55°C (-4°F to 131°F) Transmitter relative humidity: up to 95% with no condensation Transmitter pressure: 57.6 kPa to 106 kPa (8.4 psi to 15.4 psi) Charger temperature: -10°C to 50°C (14°F to 122°F) Charger relative humidity: 10% to 95% with no condensation
Battery life	Transmitter: Seven days of continuous glucose monitoring immediately following a full charge. Charger: The charger uses one new AAA battery to charge the transmitter.
Transmitter frequency	2.4 GHz band, Bluetooth® wireless technology (version 4.0)
Effective radiated power (ERP)	-12.05 dBm (0.06 mW)
Effective isotropic radiated power (EIRP)	-9.9 dBm (0.1 mW)
Operating range	Up to 1.8 meters (6 feet) in free-air
Transmitter expected service life	The transmitter expected service life is one year depending on patient usage.

Transmitter wireless communication

Quality of service

The transmitter and display device connect via a Bluetooth® low-energy technology network. The transmitter sends glucose data and system-related alerts to the display device, which verifies the integrity of received data after wireless transmission. Quality of the connection is in accordance with the Bluetooth® Specification v4.0.

Data security

The transmitter is designed to only accept radio frequency (RF) communications from recognized and linked devices. The transmitter must be paired before the display device will accept information from the transmitter.

Display devices and system components (meters and transmitters) ensure data security via proprietary means and data integrity using error checking processes, such as cyclic redundancy checks.

Traveling by air

The transmitter is safe for use on commercial airlines. Because travel rules are subject to change, it is advisable to check with the Transportation Safety Administration (TSA) before traveling.

Guidance and manufacturer's declaration

Guidance and Manufacturer's Declaration - Electromagnetic Emissions		
Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF emissions CISPR 11	CISPR 11 Group 1, Class B	The transmitter uses RF energy only for system communications. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment. Note: The preceding statement is required by IEC 60601-1-2 for Group 1, Class B devices. Since the transmitter is battery powered, its emissions will not be affected by the establishment power supply and there is no evidence of any issues associated with the use of the system in domestic establishments.
Harmonic emissions IEC 61000-3-2	Not applicable	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Not applicable	

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	For use in a typical domestic, commercial, or hospital environment.
Conducted disturbances induced by RF fields	3 V _{RMS} 150 kHz to 80 MHz 6 V _{RMS} ISM bands between 150 kHz to 80 MHz	Not applicable	Requirement does not apply to this battery powered device.
Electrical fast transient/burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	Line to Line: ±0.5 kV, ±1 kV Line to Ground: ±0.5 kV, ±1 kV, ±2 kV	Not applicable	Requirement does not apply to this battery powered device.
Note: U _T is the a.c. mains voltage prior to application of the test level.			
Voltage dips, short interruptions, and voltage variations on power supply lines IEC 61000-4-11	0% U _T ; 0.5 cycles (at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°) 0% U _T ; 1 cycle (at 0°) 70% for 25/30 cycles (at 0°) 0% for 250/300 cycles	Not applicable	Requirement does not apply to this battery powered device.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For use in a typical domestic, commercial, or hospital environment.
Proximity fields from RF wireless communications equipment IEC 61000-4-3	IEC 60601-1-2:2014, Table 9	IEC 60601-1-2:2014, Table 9	For use in a typical domestic, commercial, or hospital environment.
Note: U _T is the a.c. mains voltage prior to application of the test level.			
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 6 GHz	Portable and mobile RF communications equipment

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
	80% AM at 1 kHz	80% AM at 1 kHz	should be used no closer to any part of the transmitter than the recommended separation distance of 30 cm (12 in). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.			

Warranty

Medtronic MiniMed, Inc. (or such other legal entity as may be referred to as manufacturer on the labeling of this device "Medtronic MiniMed") warrants the Medtronic transmitter to the purchaser of the product against defects in material and workmanship for a period of one (1) year and the charger for up to one (1) year from the date of purchase.

During the warranty period, Medtronic MiniMed will replace or repair, at its discretion, any defective transmitter or charger, subject to the conditions and exclusions stated herein. This warranty applies only to new devices. In the event a transmitter or charger is replaced, the warranty period will not be extended past its original expiration date.

This warranty is valid only if the Medtronic transmitter or charger is used in accordance with the manufacturer's instructions. Without limitation, this warranty will not apply:

- If damage results from changes or modifications made to the transmitter or charger by the user, or third persons, after the date of purchase.
- If damage results from service or repairs performed by any person or entity other than the manufacturer.
- If damage results from a *Force Majeure* or other event beyond the control of the manufacturer.
- If damage results from negligence or improper use, including but not limited to: improper storage, submersion in water, physical abuse, (such as dropping).
- If damage results from use of the device in a manner other than according to the manufacturer's product labeling, instructions for use, or regulatory notifications.

This warranty shall be personal to the original purchaser. Any sale, rental or other transfer or use of the product covered by this warranty to or by a user other than the original purchaser shall cause this warranty to immediately terminate. This warranty does not apply to glucose sensors and other accessories.

The remedies provided for in this warranty are the exclusive remedies available for any breach hereof. Neither Medtronic MiniMed nor its suppliers or distributors shall be liable for any incidental, consequential, or special damage of any nature or kind caused by or arising out of a defect in the product.

All other conditions and warranties, other than mandatory statutory warranties, expressed or implied, are excluded, including the warranties of merchantability and fitness for a particular purpose.

This warranty gives the purchaser specific legal rights, and the purchaser may also have other rights that vary under local law. This warranty does not affect the purchaser's statutory rights.

Open Source Software (OSS) Disclosure

This document identifies the Open Source Software that may be separately called, executed, linked, affiliated, or otherwise utilized by this product.

Such Open Source Software is licensed to users subject to the terms and conditions of the separate software license agreement for such Open Source Software.

Use of the Open Source Software by you shall be governed entirely by the terms and conditions of such license.

The source/object code and applicable license for the Open Source Software can be obtained at the following site: <http://www.ouah.org/ogay/hmac/>.

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
Guardian 4

Introduction

The Guardian 4 transmitter (MMT-7841) with Bluetooth® wireless technology is a component of the continuous glucose monitoring (CGM) system and is compatible with the MiniMed 780G system that uses the Guardian 4 sensor (MMT-7040).

Indications for use

The Guardian 4 transmitter (MMT-7841) is a rechargeable device that powers the glucose sensor. The transmitter collects and calculates sensor data and sends the data via Bluetooth® wireless technology to the MiniMed 780G system for the management of diabetes mellitus. The transmitter is only compatible with the Guardian 4 sensor (MMT-7040) and is indicated for single-patient or multiple-patient use. The transmitter is intended for persons ages 7 years and older.

Contraindications

No contraindications are associated with Guardian 4 transmitter use. For contraindications related to CGM, see the MiniMed 780G System User Guide.

Clinical Benefits

The Guardian 4 transmitter is a component of the CGM system that provides sensor glucose values. See the MiniMed 780G System User Guide for the clinical benefits of systems that use the Guardian 4 transmitter.

User Safety

Warnings

- Always refer to the Guardian 4 Sensor User Guide for all precautions, warnings, and instructions related to the sensor. Not referring to the Guardian 4 Sensor User Guide can result in serious injury or damage to the sensor.
- Do not allow children to put small parts in their mouth. This product may pose a choking hazard that can result in serious injury or death.
- Do not change or modify the device unless expressly approved by Medtronic Diabetes. Modifying the device can cause serious injury, interfere with the ability to operate the device, and void the warranty.
- Do not expose the transmitter to Magnetic Resonance Imaging (MRI) equipment, diathermy devices, or other devices that generate strong magnetic fields (for example x-ray, CT scan or other types of radiation). Exposure to a strong magnetic field has not been evaluated and can cause the

device to malfunction, result in serious injury, or be unsafe. If the transmitter is exposed to a strong magnetic field, discontinue use and contact a local Medtronic support representative for further assistance.

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Precautions

- Do not use the transmitter adjacent to other electrical equipment that may cause interference with the normal system operation.
- Only use the Guardian 4 sensor (MMT-7040) with the transmitter. Do not use any other sensor. Other sensors are not intended for use with the transmitter and will damage the transmitter and the sensor.
- Only use the green colored tester (MMT-7736L) with the transmitter. Pockets on the transmitter are visible when connected to the tester. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.


Figure 1. Transmitter pockets

- Always use the tester when cleaning the transmitter. Do not use any other test plug with the transmitter. Use of another test plug can allow water to get into the transmitter or can prevent proper cleaning. Water can damage the transmitter.

- Do not twist the tester or sensor while attached to the transmitter. Twisting the tester or sensor will damage the transmitter.
- Do not allow the tester to come in contact with any liquid when not connected to the transmitter. A wet tester can damage the transmitter.
- Do not allow the transmitter to come in contact with any liquid when not connected to a sensor or to the tester. Moisture will damage the transmitter and a wet transmitter can damage the sensor.
- Do not clean the O-rings on the tester with any substances. Cleaning the O-rings can damage the tester.

Figure 2. O-rings

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1. Special Precautions regarding Electromagnetic Compatibility (EMC): This body worn device is intended to be operated within a reasonable residential, domestic, public or work environment where common levels of radiated "E" (V/m) or "H" fields (A/m) exist, such as cellular phones, Wi-Fi™, Bluetooth® wireless technology, electric can openers, microwave and induction ovens. This device generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the provided instructions, may cause harmful interference to radio communications.
2. Portable and mobile RF communications equipment can affect medical electrical equipment. If you encounter RF interference from a mobile or stationary RF transmitter, move away from the RF transmitter that is causing the interference.
3. Be careful when using the transmitter closer than 30 cm (12 in) to portable radio frequency (RF) equipment or electrical equipment. If the transmitter must be used next to portable RF

equipment or electrical equipment, observe the transmitter to verify correct system operation. Degradation of the performance of the transmitter could result.

4. The essential performance (EP) of the transmitter is to measure and transmit to a monitoring device the sensing device's signal value(s) within the transmitter's accuracy requirements under the specified use conditions outlined in the MiniMed 780G System User Guide and for the duration of the expected service life. If the transmitter experiences electromagnetic disturbances, either no or incorrect data may be transmitted. In such situations, refer to the operation, maintenance, and troubleshooting instructions within the applicable user guides. You may also use the tester to test if the transmitter is operating properly. If the transmitter is damaged or if it cannot communicate with the pump, contact your local Medtronic support representative for assistance.

Assistance

Contact a local Medtronic support representative for assistance and if a copy of the MiniMed 780G System User Guide is needed.

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
Pairing the transmitter

The transmitter must be paired to the system before a sensor can be used. Always refer to the MiniMed 780G System User Guide for instructions on how to pair the transmitter to the system.

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

Before proceeding, have the MiniMed 780G System User Guide available.

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
5. When the transmitter light flashes green after connecting to the sensor, use the system to start the sensor. For more instructions, see the MiniMed 780G System User Guide.
6. Attach the adhesive tab of the sensor to the transmitter.
7. After the transmitter is connected, consult the Guardian 4 Sensor User Guide for instructions on how to apply a second tape.
8. Refer to the MiniMed 780G System User Guide for pairing.

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The transmitter can be reconnected to the sensor currently in use. Simply connect the transmitter to the sensor that is already inserted. Confirm Reconnect Sensor when the pump detects the transmitter. It may take a few seconds to establish a connection when reconnecting a sensor. Reattach the adhesive tab of the sensor to the transmitter and reapply any required tape. When reconnected, the sensor goes through another warm-up period.

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

CAUTION: Only use the green colored tester (MMT-7736L) with the transmitter. Pockets on the transmitter are visible when connected to the tester. Do not use any other test plug. Other test plugs are not intended for use with the transmitter and will damage the transmitter and the tester.

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Figure 4. Transmitter components

Look inside the connector opening of the transmitter to make sure that the connector pins are not damaged or corroded. If the connector pins are damaged or corroded, the transmitter cannot communicate with the charger or pump. Contact a local Medtronic support representative. It may be time to replace your transmitter.

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
Before proceeding, have the MiniMed 780G System User Guide available.

To connect the tester:

1. Hold the transmitter and the tester as shown. Line up the flat side of the tester with the flat side of the transmitter.

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3. To test the transmitter, check the sensor icon in the app to ensure that the transmitter is sending a signal (see the MiniMed 780G System User Guide).
4. To clean the transmitter, see *Cleaning the transmitter*, page 22.
5. After testing or cleaning, disconnect the tester from the transmitter.

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

10. Place the transmitter and tester on a clean, dry cloth and air dry them completely.
11. Disconnect the tester from the transmitter by gently squeezing the arms of the tester.

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Troubleshooting

The table shown contains troubleshooting information for the transmitter, charger, and tester. For more information about troubleshooting, see the MiniMed 780G System User Guide.

Table 1. Troubleshooting issues

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Problem	Likely Cause(s)	Resolution
		tative. It may be time to replace the transmitter.
When connected to the sensor, the green light on the transmitter does not flash.	The transmitter is not fully connected. The transmitter is low on power. The sensor is not properly inserted into the body.	1. Disconnect the transmitter from the sensor. 2. Wait for five seconds and reconnect them. If the green light still does not flash, proceed to step 3. 3. Fully charge the transmitter and connect it to the tester. If the green light still does not flash, see troubleshooting on "When connected to the tester, the green light on the transmitter does not flash." If the green light flashes, proceed to step 4. 4. Disconnect the transmitter from the tester, wait at least five seconds, and connect the transmitter to the sensor. If the green light still does not flash, proceed to step 5. 5. The sensor may not be properly inserted into the body. Remove the sensor from the body and insert a new sensor.
When connected to the tester, the green light on the transmitter does not flash.	The transmitter is not fully connected. The transmitter is low on power.	1. Check the connection between the transmitter and the tester. If the green light still does not flash, proceed to step 2. 2. Fully charge the transmitter. 3. Test the transmitter with the tester again. If there is no flashing green light, call a local Medtronic support representative. It may be time to replace the transmitter.
The transmitter battery does not last for seven days.	The transmitter is not fully charged when connected to the sensor. The transmitter and pump frequently lose wireless connection.	1. Fully charge the transmitter before connecting it to the sensor. If the transmitter battery still does not last for the duration of one sensor use, proceed to step 2. 2. Move away from any equipment that can cause RF interference. For more informa-

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Problem	Likely Cause(s)	Resolution
		tion on RF interference, see the Radio Compliance Information sheet included with the pump. 3. Make sure the pump and the transmitter are located on the same side of the body to minimize any RF interference. If a fully charged transmitter battery continues to lose power before a full seven days, call a local Medtronic support representative. It may be time to replace the transmitter.
The transmitter loses connection with the pump.	The pump is out of range. There is RF interference from other devices.	1. Move away from any equipment that can cause RF interference. For more information on RF interference, see the Radio Compliance information sheet included with the transmitter. If the transmitter is still not communicating with the pump, proceed to step 2. 2. Make sure the pump and the transmitter are located on the same side of the body to minimize any RF interference. If the transmitter is still not communicating with the pump, call a local Medtronic support representative for assistance.
Note: An alarm or alert occurs and a message appears when the transmitter has lost connection with the pump for 30 minutes.		

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Quality of service

The transmitter and pump connect via a Bluetooth® low-energy technology network. The transmitter sends glucose data and system-related alerts to the pump, which verifies the integrity of received data after wireless transmission. Quality of the connection is in accordance with the Bluetooth® Specification v4.0.

Data security

The transmitter is designed to only accept radio frequency (RF) communications from recognized and linked devices. The transmitter must be paired before the pump will accept information from the transmitter.

The pump and system components (meters and transmitters) ensure data security via proprietary means and data integrity using error checking processes, such as cyclic redundancy checks.

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Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	± 8 kV contact ± 2 kV, ± 4 kV, ± 8 kV, ± 15 kV air	For use in a typical domestic, commercial, or hospital environment.
Conducted disturbances induced by RF fields	3 V _{RMS} 150 kHz to 80 MHz 6 V _{RMS} ISM bands between 150 kHz to 80 MHz	Not applicable	Requirement does not apply to this battery powered device.
Electrical fast transient/burst IEC 61000-4-4	± 2 kV 100 kHz repetition frequency	Not applicable	Requirement does not apply to this battery powered device.
Surge IEC 61000-4-5	Line to Line: ± 0.5 kV, ± 1 kV Line to Ground: ± 0.5 kV, ± 1 kV, ± 2 kV	Not applicable	Requirement does not apply to this battery powered device.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Voltage dips, short interruptions, and voltage variations on power supply lines IEC 61000-4-11	0% U_T ; 0.5 cycles (at 0°, 45°, 90°, 135°, 180°, 225°, 270°, and 315°) 0% U_T ; 1 cycle (at 0°) 70% for 25/30 cycles (at 0°) 0% for 250/300 cycles	Not applicable	Requirement does not apply to this battery powered device.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	30 A/m	30 A/m	For use in a typical domestic, commercial, or hospital environment.
Proximity fields from RF wireless communications equipment IEC 61000-4-3	IEC 60601-1-2:2014, Table 9	IEC 60601-1-2:2014, Table 9	For use in a typical domestic, commercial, or hospital environment.
Note: U_T is the a.c. mains voltage prior to application of the test level.			
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.7 GHz	10 V/m 80 MHz to 6 GHz	Portable and mobile RF communications equipment

Guidance and Manufacturer's Declaration - Electromagnetic Immunity			
Immunity Test	IEC 60601-1-2:2014 Test Level	Max foreseeable use condition per IEC 60601-1-2:2014	Electromagnetic Environment Guidance
	80% AM at 1 kHz	80% AM at 1 kHz	should be used no closer to any part of the transmitter than the recommended separation distance of 30 cm (12 in). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, should be less than the compliance level in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol:
			
Note: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.			

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
Guardian

Introduction

The Guardian application (app) is a component of the Guardian 4 system that helps manage diabetes. The Guardian 4 system has the following features:

- Records glucose values throughout the day and night
- Displays glucose values in a convenient and discreet manner on a mobile device
- Sends alerts for glucose events on a mobile device
- Shows the effects that diet, exercise, and insulin can have on glucose levels

This user guide is designed to help understand the setup and operation of the Guardian app.

Note: The user guide contains some instructions that are specific to the Android™ or iOS™ platform. Where applicable, follow the instructions specific to platform in use.

System description

The Guardian 4 system includes the following components: Guardian iOS™ app (MMT-8200) or Guardian Android™ app (MMT-8201), Guardian 4 transmitter (MMT-7841Q), Guardian 4 sensor (MMT-7040Q), and compatible tester and charger.

Continuous glucose monitoring (CGM) is a technology that lets patients continuously view glucose values. The Guardian 4 system uses a glucose sensor, the Guardian 4 sensor, placed below the skin, to continuously measure the amount of glucose in the interstitial fluid. The Guardian 4 transmitter collects and processes these glucose measurements, which are then displayed on the Guardian app. The Guardian app can also provide alerts based on sensor glucose (SG) levels.

The Guardian app also provides a user interface for entering data such as exercise, insulin, meals, blood glucose (BG) values, and uploading information to the CareLink Personal website.

The Guardian app is available to download from the Apple™ App Store™ or the Google Play™ store and requires the Guardian 4 sensor and Guardian 4 transmitter to function.

Note: This product should only be used with supported mobile devices and operating systems. Refer to the local Medtronic website for information about supported devices and operating systems.

Diabetes treatment decisions

Use a BG meter reading to make treatment decisions if no sensor data is available or if unsure that SG values are correct.

Indications for use

The Guardian app is intended to be used with the Guardian 4 system.

The Guardian 4 system is a real-time, CGM system indicated for the management of diabetes in persons age 7 years and older.

The Guardian 4 system does not require calibration and is designed to replace finger stick BG testing for diabetes treatment decisions. Treatment decisions should be made based on the combination of SG readings and trend arrows. The system can track trends in glucose concentrations, and aid in detecting events of hyperglycemia and hypoglycemia and helping both current and long-term therapy decisions.

The system is intended for use only by patients and caregivers using a compatible mobile device and operating system, and who have sufficient experience to adjust mobile device audio and notification settings.

Contraindications

Successful CGM use requires sufficient vision or hearing ability to allow recognition of the alerts generated by the Guardian app.

Clinical benefits

The Guardian app is software that provides information that is used for diabetes management but does not provide any direct therapy. As a result, the app does not have direct therapeutic benefits but can support the management of diabetes by patients and caregivers.

User safety

Warnings

Follow the instructions and safety warnings in this user guide to receive alerts as intended. Missing alerts from the Guardian app may result in undetected low and high glucose levels.

- CGM is not recommended for people who are unwilling or unable to perform BG tests as required or for people who are unwilling or unable to maintain contact with their healthcare professional.
- Always allow notifications for the Guardian app. If notifications are turned off, the app will not send any alerts, including Urgent Low Alert.

- Do not use the Guardian app without understanding how the mobile device settings work. If the mobile device settings are not set up correctly, the app may not send SG alerts.
- Make sure Bluetooth™ is on, even if the mobile device is in airplane mode. If Bluetooth™ is off, the app will not send SG information or alerts.
- Do not use the Guardian app if the mobile device screen or speakers are damaged. If the mobile device is damaged, the app may not send SG alerts and SG information may not be shown correctly.
- Do not force close the Guardian app. If the app is closed, the app will not send SG information or alerts.
- Check the Guardian app occasionally to make sure it is running. The mobile device may close the Guardian app automatically when another app is in use, such as a game. If the Guardian app is closed, the app will not send SG alerts.
- Do not let the mobile device shut down due to low battery, or the app will not send SG alerts. Use of the app may deplete the mobile device battery more quickly. Have a charger available to charge the battery if needed.
- Always make sure to open the app after the mobile device restarts to ensure the app sends SG alerts.
- Make sure to set the snooze to a short enough time so that the app sends an alert again if glucose levels do not improve. When an SG alert is snoozed, the app will not send that alert again during the length of the set snooze time.
- Do not use continuous glucose monitoring if hydroxyurea, also known as hydroxycarbamide, is taken. Hydroxyurea is used to treat certain diseases, such as cancer and sickle cell anemia. Hydroxyurea use results in higher sensor glucose readings compared to blood glucose readings. Taking hydroxyurea while using continuous glucose monitoring can result in inaccurate or missed alerts, and substantially higher sensor glucose readings in reports than actual blood glucose readings.
- Always check the label of any medication being taken to confirm if hydroxyurea or hydroxycarbamide is an active ingredient. If hydroxyurea is taken, consult a healthcare professional. Use additional blood glucose meter readings to verify glucose levels.
- Taking medications that contain paracetamol or acetaminophen, such as fever reducers or cold medicine, while wearing the sensor may falsely raise SG readings. The level and duration of sensor inaccuracy depends on how much medication is working in the body and will be different for each person. If medication is taken with paracetamol or acetaminophen, use a BG meter reading for dosing decisions in order to avoid a possible hypoglycemic event. Consult a healthcare professional when determining insulin dosage.
- Turning off automatic updates on the mobile device may help to avoid unintentionally updating to an operating system that is not confirmed as compatible with the Guardian app.

- Do not root or jailbreak the mobile device. Rooting the Android™ device or jailbreaking the iOS™ device means to change the software in a way the manufacturer did not intend. If the mobile device is changed in this way, the Guardian app will display an error message when launched and will not continue to operate.

Android™ Users

- Allow **Do Not Disturb Permissions** and **Notifications** for the Guardian app. If **Do Not Disturb Permissions** or **Notifications** are turned off, the app will not send any alerts, including Urgent Low Alert.
- The Android™ device may include battery saver settings. If these settings are active for the Guardian app, the mobile device will shut down the app and prevent it from sending alerts and notifications. Check the mobile device periodically to confirm that the Guardian app is active.
- The **Digital Wellbeing** feature is intended to prevent alerts. If this feature is enabled, the Guardian app will not send alerts. If the app timer is set in the **Digital Wellbeing** feature for the Guardian app, the app will shut down and will not send any SG alerts when the timer expires.
- Do not lower the vibration level to the lowest setting in the Android™ **Vibration intensity** menu. If vibration is lowered to the lowest setting, the app will not send any vibrations with the alerts, including Urgent Low Alert.
- Alerts for the Guardian app will sound through headphones when headphones are connected. If headphones are connected when not in use, SG alerts may not be heard.

iOS™ Users

- Allow **Critical Alerts** and **Notifications** for the Guardian app. If **Critical Alerts** or **Notifications** is turned off, the app will not send any alerts, including Urgent Low Alert.
- The **Downtime** setting within the **Screen Time** feature is intended to prevent alerts during the **Downtime** period. If this setting is enabled, the Guardian app will be prevented from sending alerts.
- Do not turn off vibration in the iOS™ Accessibility menu. If vibration is turned off, the app will not send any vibrations with the alerts, including Urgent Low Alert.
- If a serious incident related to the device occurs, immediately report the incident to the manufacturer and local competent authority.

Precautions

Diabetes treatment decisions should be made based on a combination of SG readings, trend arrows, glucose target ranges, active alerts, and recent events (such as insulin doses, exercise, meals, and medications).

Using CGM information to make treatment decisions

After becoming familiar with CGM, treatment decisions should be made based on all the information available, including the following:

- SG readings
- Trend arrows
- Active SG alerts
- Recent events such as insulin doses, medication, meals, exercise, etc.

Refer to a healthcare professional for glucose target ranges.

Risks and side effects

Software malfunction or failure to follow instructions of the Guardian app can lead to hypoglycemia, hyperglycemia, diabetic ketoacidosis, or possibly death.

Assistance

Contact the local representative or refer to the local Medtronic website for assistance. Refer to the Medtronic Diabetes International Contacts list in this user guide for contact information.

Guardian app setup

To set up the Guardian app, follow the instructions on the screen.

Pairing the transmitter

The transmitter must be paired to the system before a sensor can be used. Follow the instructions on the screen to pair the transmitter to the mobile device.

New sensor setup

Follow the video instructions on how to insert the sensor. Then follow the on-screen instructions to complete the sensor start up.

Connecting the transmitter to the sensor


Refer to the *Guardian 4 Transmitter User Guide* for instructions on how to connect the transmitter to the sensor.

Completing the app setup

Continue to follow the on-screen instructions to enable notifications. For more information, see *Alert settings, page 19*.

Home screen

The following figure shows the Home screen of the app.

Note: This screen may vary depending on the mobile device and platform.


Note: The Home screen does not display sensor information the first time the app is opened. The first SG reading appears after the transmitter is successfully paired and the sensor warm-up is complete.

Item	Description
Trend arrows	Display the glucose trend and the rate at which the most recent SG level has risen or fallen. For details on how to set the rate of change for falling and rising SG values, see <i>Setting rate alerts</i> , page 22.
Menu screen	Provides access to the menu screen, which displays the Guardian 4 system status information and allows the user to define alert settings, set up the sensor, calibrate, pair the transmitter, access CareLink Personal settings, view the logbook, and access help documentation. For details, see <i>Menu screen</i> , page 15.
Continuous SG trace	Displays current and previous SG readings. Tap any point on the graph to view the details of the selected SG reading or event. For more navigating tips, see <i>Navigating through the SG graph</i> , page 24.
System status icons	Display the corresponding status icon if the transmitter battery is empty, the sensor life has expired, the connection with the transmitter is lost, notifications are silenced, or the status of these features is unknown. For details, see <i>System status icons</i> , page 16.
Time interval	Displays preset time intervals of 3, 6, 12, and 24 hours. Tap this icon to switch between time intervals.
Current SG reading	Displays the current SG reading, which the transmitter calculates and sends wirelessly to the app.
High and low SG limits	Display a line showing the high and low SG limits on the sensor graph. The orange line indicates the high SG limit; the red line indicates the low SG limit.
Sensor status message	Displays the latest active status notification. If a current SG reading is not available, the sensor status message will appear where the SG information would be displayed.
SG graph	Displays a white line representing SG readings over a selected time interval. It also indicates high and low SG limits. For details on the sensor graph, see <i>Sensor graph</i> , page 24.

Item	Description
Day indicator	Displays the date of the data shown in the graph. Scroll further back into historical data to see Today, Yesterday, and days of the previous week.
Add event	Displays the events screen to enter additional information, such as exercise, BG readings, meals, or insulin taken. Certain BG readings entered here may be also used for sensor calibration. For details, see <i>Events, page 25</i> .

Menu screen

Tap  in the upper left hand corner of the Home screen to see the menu screen.

Item	Description
Close (X)	Closes the menu screen and returns to the Home screen.
System status bar	Displays icons that show a quick status of the Guardian 4 system, including sensor life, transmitter battery, and transmitter communication. For details on the individual icons and their status, see <i>System status icons, page 16</i> .
Alert Volume	Allows the option to adjust the volume of alerts. Alerts will sound louder and repeat if not acknowledged. For details, see <i>Alert settings, page 19</i> .
Alert Volume Information	Displays a description of Alert Volume slider and Mute All Alerts button. For details, see <i>Alert settings, page 19</i> .
Mute All Alerts	Allows the option to mute all alerts for a specified duration. Urgent Low Alert will always vibrate. For details, see <i>Alert settings, page 19</i> . Note: Urgent Low Alert will not vibrate if vibration is turned off in the iOS™ Accessibility menu or lowered to the lowest setting in the Android™ Vibration intensity menu.
Alert Settings	Display the options to set glucose monitoring alerts. For details, see <i>Alert settings, page 19</i> .
Sensor & Transmitter	Displays the instructions to start a new sensor, calibrate, and pair a transmitter. For details, see <i>New sensor setup, page 12</i> .
CareLink	Provides access to a menu to change the CareLink user login information. For details, see <i>Syncing data to the CareLink Personal website, page 32</i> .
Logbook	Displays a history of alerts and events that have occurred in the Guardian 4 system. For details, see <i>Logbook, page 29</i> .
Help	Provides the user with software version information, access to the app user guide, getting started guide, setup wizard and the end user license agreement.

System status icons

The system status icons appear at the top of the menu screen. These icons provide a way to quickly check the status of the system. If any condition becomes critical and requires immediate attention, the corresponding status icon is also displayed on the Guardian Home screen. The icons displayed on the Home screen are interactive and provide more information about the current status.

Icon name	Description
Transmitter battery	As the battery life is used, the icon changes in the following order: 

Icon name	Description
	indicates that approximately 100% of the battery capacity remains, which means at least 7 days of use remain indicates that approximately 50% of the battery capacity remains indicates that up to one day of use remains indicates that the battery is empty indicates that the battery status is unknown
Sensor life	As the sensor life is used, the number on the icon indicates the number of days left before the sensor needs to be replaced, and the color of the icon changes as shown: When a new sensor is inserted, the icon is solid green. When there is less than one day left, the 1-day icon remains red, and the text below it indicates the number of hours that remain. The question mark indicates that the sensor life is unknown. When the sensor is expired, an "X" is displayed.
Transmitter communication	Displays the communication status between the transmitter and the app: indicates the transmitter is active and connected indicates there is a communication error, the transmitter is not paired to the mobile device, or Bluetooth™* is off

Alerts

The Guardian app provides **glucose alerts** and **system status alerts**. These alerts provide information about glucose levels and the status of the Guardian 4 system.

Glucose alerts

Set **glucose alerts** for the following situations:

- SG values are rising or falling faster than the set limits
- SG values have gone above or below the set limits
- SG values are predicted to go above or below the set limits

The **glucose alerts** in the Guardian 4 system are listed in the following table:

Alert type	Description
High Alert	SG values have gone above the set high limit.
High predicted	The SG value is predicted to go above the set high limit, within a set period of time (up to 60 minutes ahead).
Rise Alert	The SG value is rising faster than the set rate (corresponding to the rising arrows displayed next to the SG level).
Low Alert	The SG value has gone below the set low limit.
Low predicted	The SG value is predicted to go below the set low limit, within a set period of time (up to 60 minutes ahead).
Fall Alert	The SG value is falling faster than the set rate (corresponding to the falling arrows displayed next to the SG level).
Urgent Low Alert	The SG level has gone below 3.0 mmol/L (54 mg/dL).

Note: If the Guardian app is in "Alerts Disabled", the app will not send any alerts, including the Urgent Low Alert.

Glucose alerts can be customized and are set up as described in *Alert settings, page 19*.

System status alerts

The Guardian 4 system also has system **status alerts** that provide information about actions needed to ensure the correct functioning of the system. See the status alerts table below for a complete listing of these alerts. For more information on how to address these alerts, see *Troubleshooting, page 33*.

The **status alerts** in the Guardian 4 system are listed below:

Make sure to act on these alerts to continue to receive SG information in the future.

Alert type	Description
Mobile device battery low	The mobile device battery has reached or fallen below 20% of its power.
Calibrate now	Measure the BG value with the meter and calibrate the sensor to get SG readings.
Calibration not accepted	The BG meter value could not be used to calibrate.
Change sensor	The sensor is not working properly and needs to be replaced.
Lost communication	The Guardian app and transmitter have not been communicating for 30 minutes. The app may have closed if there are too many apps running at the same time or if there is radio frequency interference.

Alert type	Description
Sensor end of life	Sensor has reached its maximum life of 7 full days.
Sensor updating	The sensor is updating and no SG information will be available. Measure the BG value with the meter.
Transmitter battery empty	The transmitter battery is empty and needs to be recharged. The app is no longer receiving sensor information.
Transmitter error	The transmitter is trying to fix a problem.
Jailbroken device detected	The software on the iOS™ mobile device has been changed such that it no longer works in the way the manufacturer intended. The app cannot be used with a jailbroken device. Changing the mobile device operating software causes the app to stop working.
Rooted device detected	The software on the Android™ mobile device has been changed such that it no longer works in the way the manufacturer intended. The app cannot be used with a rooted device. Changing the mobile device operating software causes the app to stop working.

Alert settings

Setting low alerts

Low alert settings include the following:

Low alert setting	Description
Day Starts at	Set when the day begins, what the low limit is during the day, and when the app should send alerts at this setting.
Night Starts at	Set when the night begins, what the low limit is during the night, and when the app should send alerts at this setting.
Low Limit	The Low Limit is the SG value on which the other low settings are based. The low limit can be set from 3.3 to 5.0 mmol/L (60 to 90 mg/dL). On the SG graph, the low limit appears as a red horizontal line at the set value.
Don't alert me	When Don't alert me is selected, the app will not send any alerts for low glucose. The app still sends an Urgent Low Alert since this alert is always on.
At Low Limit	When At Low Limit is selected, the system displays a Low Sensor Glucose Alert when the SG value reaches or falls below the low limit.
Before Low Limit	When Before Low Limit is selected, the app sends a Low predicted alert any time the SG is predicted to reach the low limit. This alert notifies of potential low glucose levels before they occur.

Low alert setting	Description
Before and at Low Limit	When Before and at Low Limit is selected, the app sends a Low predicted alert any time the SG is predicted to reach the low limit and when the SG value reaches or falls below the low limit.
Time Before Low	Time Before Low is available only when Before Low Limit or Before and at Low Limit is selected. Time Before Low determines when the app will send a Low predicted alert, if the SG values continue to decrease at the current rate of change. The option is available to set a time between ten minutes and one hour.
Max Volume at Night	Max Volume at Night is available only in the Night Time section and when Alert Me is selected. Toggle Max Volume at Night on for low alerts to make all low alerts sound at max volume at night.

To set low alerts:

1. On the Home screen, tap  and select **Alert Settings > Low Alerts**.
2. Tap **Day starts at** and set the desired start time. Tap **Save**.
3. Tap **Low Limit** and set the desired limit between 3.3 to 5.0 mmol/L (60 to 90 mg/dL). Tap **Save**.
4. To receive alerts when the SG value approaches the low limit, tap **Alert Me** and choose one of the following:
 - a. At Low Limit
 - b. Before Low Limit
 - c. Before and at Low Limit
5. If **Before Low Limit** or **Before and at Low Limit** is selected, the **Time Before Low** appears. Set the desired duration of time to receive a Low predicted alert and tap **Save**.
6. Return to the Low Alerts screen and tap **Night starts at**. Set the desired start time and tap **Save**. Repeat steps 3 through 5 to set night time alerts.
7. Toggle **Max Volume At Night** on or off. When **Max Volume At Night** is toggled on, the low alerts sound at max volume at night.
8. When the alerts are set up, tap  to return to Alert Settings screen.
9. Tap  to return to the menu screen.
10. Tap  to close the menu and return to the Home screen.

Setting high alerts

High alert settings include the following:

High alert setting	Description
Day Starts at	Set when the day begins, what the high limit is during the day, and when the app should send alerts at this setting.
Night Starts at	Set when the night begins, what the high limit is during the night, and when the app should send alerts at this setting.
High Limit	The high limit is the SG value on which the other high settings are based. The high limit can be set from 5.5 to 22.2 mmol/L (100 to 400 mg/dL). On the SG graph, the high limit appears as an orange horizontal line at the set value.
Don't Alert Me	When Don't Alert Me is selected, the app does not send any alerts for high glucose.
At High Limit	When At High Limit is selected, the system displays a High Sensor Glucose alert when the SG value reaches or exceeds the high limit.
Before High Limit	When Before High Limit is selected, the app sends a High predicted alert any time the SG value is predicted to reach the high limit. This notifies of potential high glucose levels before they occur.
Before and at High Limit	When Before and at High Limit is selected, the app sends a High predicted alert any time the SG is predicted to reach the high limit and when the SG value reaches or rises above the high limit.
Time Before High	The Time Before High option is available only when Before High Limit or Before and At High Limit is selected. This option determines when the app sends a High predicted alert if the SG values continue to increase at the current rate of change. Set a time between ten minutes and one hour.
Max Volume at Night	Max Volume at Night is available only in the Night time section and when Alert me is selected. Toggle Max Volume at Night on for high alerts to make all high alerts sound at max volume at night.

To set high alerts:

1. On the Home screen, tap  and select **Alert Settings > High Alerts**.
2. Tap **Day starts at** and set the desired start time. Tap **Save**.
3. Tap **High Limit** and set the desired limit between 5.5 to 22.2 mmol/L (100 to 400 mg/dL). Tap **Save**.
4. To receive alerts when the SG value approaches the high limit, tap **Alert Me** and choose one of the following:
 - a. At High Limit
 - b. Before High Limit
 - c. Before and at High Limit

5. If **Before High Limit** or **Before and at High Limit** is selected, the **Time Before High** appears. Set the desired time to receive a High predicted alert and tap **Save**.
6. Return to the High Alerts screen and tap **Night starts at**. Set the desired start time and tap **Save**. Repeat steps 3 through 5 to set night time alerts.
7. Toggle **Max Volume At Night** on or off. When **Max Volume At Night** is toggled on, high alerts sound at max volume at night.
8. When the alerts are set up, tap  to return to Alert Settings screen.
9. Tap  to return to the Menu screen.
10. Tap  to close the menu and return to the Home screen.

Setting rate alerts

Rate alerts notify when the SG value is rising (Rise Alert) or falling (Fall Alert) equal to or faster than the set rate. These alerts help to understand how glucose levels are affected, for example, by meals or exercise.

On the Home screen, these fall or rise rates are indicated by arrows, as shown in *Home screen, page 13*. The more arrows, the faster the rate of change.

↑	SG has been rising at a rate of 0.06 mmol/L (1 mg/dL) or more per minute, but less than 0.11 mmol/L (2 mg/dL) per minute.
↓	SG has been falling at a rate of 0.06 mmol/L (1 mg/dL) or more per minute, but less than 0.11 mmol/L (2 mg/dL) per minute.
↑↑	SG has been rising at a rate of 0.11 mmol/L (2 mg/dL) or more per minute, but less than 0.17 mmol/L (3 mg/dL) per minute.
↓↓	SG has been falling at a rate of 0.11 mmol/L (2 mg/dL) or more per minute, but less than 0.17 mmol/L (3 mg/dL) per minute.
↑↑↑	SG has been rising at a rate of 0.17 mmol/L (3 mg/dL) or more per minute.
↓↓↓	SG has been falling at a rate of 0.17 mmol/L (3 mg/dL) or more per minute.

To set rate alerts:

1. On the Home screen, tap  and select **Alert Settings > Fall & Rise Alerts**.
2. Toggle **Fall Alert** on.
3. Tap the arrow option with the desired fall rate.
4. Toggle **Rise Alert** on.
5. Tap the arrow option with the desired rise rate.
6. When fall and rise alerts are set up, tap  to return to Alert Settings screen.

Setting alert snooze time

The snooze feature gives the option to set a snooze time for alerts. If the alert condition still persists after the set snooze time, the app sends a notification of the alert. The set snooze time for high and rise alerts may be different from the set snooze time for low and fall alerts.

To set the alert snooze time:

1. On the Home screen, tap , and select **Alert Settings > Snooze Time**.
2. Tap **Low & Fall Alerts** and set the desired snooze time for Low & Fall Alerts. Tap **Save**.
3. Tap **High & Rise Alerts** and set the desired snooze time for High & Rise Alerts. Tap **Save**.

Adjusting alert volume

All Guardian alerts sound at the volume set in the app. Alerts sound louder and repeat if not acknowledged.

To adjust alert volume:

1. On the Home screen, tap .
2. At the top, there is a slider. Swipe left to turn the volume down. Swipe right to turn the volume up.

Note: If the volume is set at 0%, a popup message appears stating: "Alert volume set at 0%. If you don't respond to the initial alert, it will sound and repeat. Your alerts are NOT muted. To mute alerts, tap **Mute All Alerts**."

Muting alerts

All Guardian alerts will be muted for the selected time. Urgent Low Alert will always vibrate.

To mute alerts:

1. On the Home screen, tap .
2. Tap **Mute All Alerts** and 4 options appear.
 - a. 30 min
 - b. 1 hour
 - c. 4 hours (max)
 - d. Custom

If Custom is selected, adjust the mute time to the desired length and tap **Mute**.

3. A popup message appears to confirm the alerts are muted. To cancel the mute, tap **Cancel Mute**.

Acting on Guardian app alerts

The Guardian app sends alerts, similar to notifications sent from other apps on the mobile device. All alerts vibrate according to the mobile device vibration settings.

Open the Guardian app to address an alert. Dismissing the notification only removes it from the list of notifications on the mobile device. If the notification is dismissed from the mobile device, but isn't acted on it in the Guardian app, the alert may repeat.

The alert is displayed on the screen when the app is opened. Click **OK** to clear system status alerts. For SG alerts, either swipe the alert up, which will snooze the alert for the set snooze period, or swipe the alert down to set a new snooze period.

Sensor graph


The sensor graph displays the current SG reading. It also gives the option to view a history of SG readings and entered events.

Navigating through the SG graph

- Swipe the center of the graph right and left to view historical data and to return to the current graph location.
- Pinch and stretch the center of the graph to zoom in and out on the graph data.
- Tap the graph twice to view the graph at the selected 3-hour, 6-hour, 12-hour, and 24-hour zoom levels.
- Tap the horizontal time axis above the graph to center the selected time point on the graph. This gives the option to view the details of the selected SG reading or event in an information box that appears above the graph. For details, see *Graph information boxes, page 24*.
- Double tap the vertical SG values axis to the right of the graph to return to the current SG reading displayed on the graph.

Graph information boxes

Tap on the graph to view information about an event or SG reading. The tapped time is marked by a vertical cursor on the graph, and a box with the information about the tapped event or SG reading appears above the graph. The following figure provides an example:

For SG readings, information boxes display the SG value, its date and time, and trend arrows, similar to the current SG value information displayed on the Home screen. If there is no SG value for the selected time point, the information box displays the sensor status message for that point.

For event markers, information boxes display the details specific to each event type. For details, see *Events, page 25*.

Drag the cursor across the graph, to skip at five-minute intervals, to pinpoint a specific SG value on the graph.

Information boxes appear for a few seconds and then close.

Events

Events help capture information that may affect glucose levels. Use the Events screen on the app to enter and save certain types of events.

Event icon	Event name	Description
	BG	The BG marker represents the BG meter readings. These can be used either to calibrate the system or simply to manage diabetes without calibrating the system.
	Insulin	The Insulin marker represents the type and amount of insulin delivered.
	Meal	The Meal marker represents the amount of carbohydrates consumed, food or drink.
	Exercise	The Exercise marker represents the intensity and duration of exercise routine.
	Notes	The Notes marker represents other information relevant to diabetes management. For example, a record of other medications taken, illness, or stress.

Note: Make it a practice to enter events when they happen. Events can be deleted and entered again.

Entering BG meter readings

The app gives the option to enter BG meter readings. For example, if measuring BG when eating or when the SG value is rising or falling rapidly, enter the BG meter reading into the app.

The app also gives the option to use the entered BG meter reading for calibration purposes, if calibration is allowed at the time when the event is entered. To use the BG meter reading for calibration, the value must be between 2.8 mmol/L and 22.2 mmol/L (50 mg/dL and 400 mg/dL.) When calibrating the sensor from the Blood Glucose screen, enter the current BG value, or a value that is up to ten minutes old.

Note: Enter a BG value between 1.1 mmol/L and 33.3 mmol/L (20 mg/dL and 600 mg/dL) into the Events screen. However, for sensor calibration, a BG value between 2.8 mmol/L and 22.2 mmol/L (50 mg/dL and 400 mg/dL) must be entered.

To enter a BG meter reading into the Events screen:

1. Measure the BG with a BG meter.
2. Tap  at the bottom of the Home screen.
3. Tap  on the Events screen.
4. To change the date or time for the entry, tap **Time** and make the appropriate changes.
5. Use the number pad to enter the BG meter reading.
6. Make sure the value displayed above the number pad is correct. If the value is incorrect, tap  to clear it and enter the correct value.
7. To avoid calibrating the sensor with this reading, tap **Save** at the bottom of the screen.
8. To calibrate the sensor with the entered BG meter reading when calibration is available, tap **Calibrate** at the bottom of the screen.

The app returns to the Home screen after selecting **Save** or **Calibrate**, and  appears on the graph at the selected time.

The following icons may appear during calibration. Select the icon to see more information about the calibration status.

 indicates that the calibration status is unknown

 indicates that calibration is currently not permitted

Entering insulin injection information

If delivering insulin using an insulin pump, insulin pen, or a syringe, use the app to record the amount of insulin injected.

To enter the type and amount of insulin injected:

1. Make a note of the amount and type of injection taken.
2. At the bottom of the Home screen, tap .
3. On the Events screen, tap .
4. To change the date or time for the entry, tap **Time** and make the appropriate changes.
5. Use the number pad to enter the insulin amount.
6. Make sure the value displayed above the number pad is correct. If not, tap  to clear it and enter the correct value.
7. Tap **Type** and select the type of insulin taken:
 - Rapid-acting
 - Long-acting
 - Unspecified
8. Tap  to return to the Insulin Event screen.

9. Tap **Save** at the bottom of the screen. The app returns to the Home screen, and ■ appears on the graph at the selected time.

Entering meal information

Use the app to record information about the carbohydrates consumed with meals or snacks.

To enter meal information:

1. Determine the total amount, in grams or exchanges, of carbohydrates in the planned meal, snack, or drink.
2. At the bottom of the Home screen, tap ●.
3. On the Events screen, tap ■.
4. To change the date or time for the entry, tap **Time** and make the appropriate changes.
5. Use the number pad to enter the amount of carbohydrates consumed.
6. Make sure the value displayed above the number pad is correct. If not, tap ■ to clear it and enter the correct value.
7. Tap **Save** at the bottom.

The app returns to the Home screen, and ■ appears on the graph at the selected time.

Entering exercise information

Use the app to enter information about exercise regimen. Make sure to be consistent and enter the marker either before or after each period of exercise.

To enter exercise information:

1. Make a note about the length of time exercised (duration) and how difficult or easy the exercise was (intensity).
2. At the bottom of the Home screen, tap ●.
3. On the Events screen, tap ■.
4. To change the date or time, tap **Time** and set the desired day and time for this exercise entry.
5. Tap **Intensity** and select **Low**, **Medium**, or **High** to indicate the intensity of exercise.
6. Tap **Duration**.
7. Tap **Save** at the bottom.

The app returns to the Home screen, and ■ appears on the graph at the selected time.

Entering notes

Use the app to enter events other than BG measurements, insulin injections, carbohydrates consumed, and exercise information. For example, record information such as when medications are taken, illness, or stress.

To enter notes:

1. At the bottom of the Home screen, tap .
2. On the Events screen, tap .
3. Use the text field to enter the relevant information.
4. Tap **Save** at the bottom.

The app returns to the Home screen, and  appears on the graph at the selected time.

Logbook

The Logbook screen displays a history of alerts and events that occurred on the selected day, with the most recent entries at the top of the list.

To view Logbook entries:

1. On the Home screen, tap  and select **Logbook**.
2. Do any of the following to view the desired information:
 - Tap **Alerts** or **Events** to filter the list by the specific type. To view the entire list, select **All**.
 - Swipe down and up on the list to view the entries.
 - To delete an event entry, swipe it left and tap **Delete**.

Note: Alerts or calibration events cannot be deleted.

- Tap the desired entry to expand it and view the details. Use  and  at the top of the expanded view to scroll through the list entries.

Before using SG readings to make treatment decisions

Before using SG readings to make treatment decisions, consult with the healthcare professional to do the following:

- Develop a diabetes management plan
- Determine personal glucose target ranges

Until becoming comfortable using SG readings to make treatment decisions, continue to use BG meter readings to make treatment decisions and confirm SG readings.

If the SG readings do not match symptoms, use a BG meter reading to confirm the SG value. If SG readings continue to be different than symptoms, consult the healthcare professional about how to use SG readings to help manage diabetes.

When to use BG meter readings

In the following conditions, use the BG meter readings to make treatment decisions.

- **The medication taken contains paracetamol or acetaminophen.**

Wait to use the SG reading to make treatment decisions until no longer taking medication. Medication that contains paracetamol or acetaminophen used to reduce fevers or treat cold and flu symptoms, can falsely raise SG readings.

- **The most recent SG reading is unavailable.**

If a new sensor is inserted or the Sensor updating message is displayed the SG reading will be unavailable. Check the BG value with a BG meter reading, and use the BG reading to make treatment decisions until sensor readings are available.

- **Symptoms do not match the SG value.**

Check the BG value with a BG meter reading before using the SG value to make treatment decisions.

Enter the BG meter reading into the app immediately after checking BG. Review the BG meter instructions for guidance on how to check the BG.

CAUTION: Only use blood from fingertips to check BG levels. Do not use any other part of the body to obtain blood to test BG.

- Always make sure hands are clean before testing BG.
- Avoid use of an old BG reading or reuse of BG readings from previous calibrations.
- If for some reason calibration is unsuccessful, wait at least 15 minutes before attempting another calibration.

Note: The Guardian app contains a feature which provides alerts prior to SG levels reaching a set high or low limit. SG readings can vary from actual BG readings, therefore there may be situations when alerts are displayed and BG levels have not reached the set high or low limit. Consult the healthcare professional for actions and adjustments to alerts. There may also be situations when BG levels have reached the set high or low limit without any alerts displayed on the Guardian app. If there are symptoms of high or low BG levels, check the BG value immediately. Refer to the healthcare professional for actions and adjustments for high or low BG levels.

Using SG readings to make treatment decisions

When using CGM, there are several things to consider to help make treatment decisions.

Look at the most recent SG value along with the SG graph, trend arrows, and alerts. The SG graph helps to understand how SG values may have recently changed. The trend arrows indicate what the SG level may be soon.

- If the SG reading is lower than the SG target, and CGM shows **↑↑**, this indicates the SG value is increasing. As a result, consider waiting to treat or adjust treatment for the low SG value. Consider symptoms before making treatment decisions based on SG value.
- If the SG reading is higher than the SG target, and CGM shows **↓↓**, this indicates the SG value is decreasing. As a result, consider waiting to treat or adjust treatment for the high SG value. Consider symptoms before making treatment decisions based on the SG value.

Before taking an insulin dose based on the current SG, consider whether insulin from a previous dose is still lowering glucose levels. Stacking insulin is the process of taking an additional dose of insulin while there is still active insulin. Insulin stacking may cause low BG levels.

The following table can help to make treatment decisions.

Arrows Displayed	Low glucose	High glucose	Target glucose
None	May need to take a fast-acting sugar.	May adjust insulin to correct a high glucose to reach the glucose target range. Do not stack insulin.	Consider last insulin dose and meal taken. Watch for any changes to glucose levels. Do not stack insulin.
1 UP Arrow ↑	Watch for any changes to glucose levels, consider symptoms, and wait to see if SG returns to target range. Make sure to not overtrear for a low SG level.	May adjust insulin to correct a high glucose to reach the glucose target range. Do not stack insulin.	Consider last insulin dose and meal taken. May need to take insulin to stay within the glucose target range. Do not stack insulin.
2 UP Arrows ↑↑	Watch for any changes to glucose levels, consider symptoms, and wait to see if	May adjust insulin to correct a high glucose to reach the glucose target range. Do not stack insulin.	May need to take insulin to stay within the glucose target range if insulin was not already taken with a recent meal or snack.

Arrows Displayed	Low glucose	High glucose	Target glucose
	SG returns to target range. Make sure to not overtreat for a low SG level.		Do not stack insulin.
3 UP arrows ↑↑↑	Watch for any changes to glucose levels, consider symptoms, and wait to see if SG returns to target range. Make sure to not overtreat for a low SG level.	May adjust insulin to correct a high glucose to reach the glucose target range. Do not stack insulin.	May need to take insulin to stay within target range. Do not stack insulin.
1 DOWN arrow ↓	May need to take a fast-acting sugar or eat a snack.	Consider last insulin dose and activity taken. May need to watch and wait to reach target range. Do not stack insulin.	May need to take a fast-acting sugar and eat a snack.
2 DOWN arrows ↓↓	May need to take a fast-acting sugar.	Consider last insulin dose and activity taken. May need to watch and wait to reach target range. Do not stack insulin.	May need to take a fast-acting sugar and eat a snack.
3 DOWN arrows ↓↓↓	May need to take a fast-acting sugar.	Consider last insulin dose and activity taken. May need to watch and wait to reach target range. Do not stack insulin.	May need to take a fast-acting sugar and eat a snack.

Syncing data to the CareLink Personal website

Guardian 4 offers remote monitoring by Care Partners and daily uploads to the CareLink Personal website. This is done using an automatic feature called Sync to CareLink.

This automatic Sync to CareLink feature sends data displayed in the app to the CareLink Connect tab of the CareLink Personal website. This data is sent approximately every five minutes when an Internet connection to the website is available. This feature also automatically sends sensor history information used to create CareLink Personal reports roughly every 24 hours. This information can be viewed by a patient or a care partner on the CareLink Personal website at carelink.minimed.eu. Please note that the

mobile device must be connected to the Internet to send data to the website. If using a cellular connection, the mobile device provider's data rates may apply.

If the Sync to CareLink toggle is turned off, the app no longer sends sensor information to the CareLink Personal website.

The **Upload Now** button gives the option to immediately send sensor history data to the website for generating reports.

Sharing CareLink Personal data with care partners

Now that the app is synced with the CareLink Personal website, users can share their data with a care partner, such as a family member or friend. Care partners can visit carelink.minimed.eu on a personal computer to create an account and send a follow request.

For more information, refer to the *CareLink Personal User Guide* that can be found on the CareLink Personal website.

Disconnecting the transmitter from the sensor

Always refer to the *Guardian 4 Transmitter User Guide* for instructions on how to disconnect the transmitter from the sensor.

Reconnecting the existing sensor

If necessary, a transmitter can be reconnected to an existing sensor. When the app detects the connection, confirm that the sensor is an existing sensor. It may take a few seconds to establish a connection when connecting an existing sensor. If an existing sensor is reconnected, the sensor will go through another warm-up period.

Troubleshooting

The following table contains troubleshooting information for the alerts.

Alerts

Problem	Likely Cause(s)	Resolution
Lost communication alert	Another app is in use, such as a game, that takes up a lot of the memory on the mobile device. This means that the Guardian app stops running and can't communicate with the transmitter.	Open the app to ensure it is running properly. Check periodically to see if the app is still running in the background in order to receive alerts and SG values.
	The app has been closed. The Guardian app has stopped running and can't communicate with the transmitter. Potential causes include using other apps and features (for example task manager apps), or selecting Force Stop for Guardian from an Android™* device Settings menu.	
Lost communication alert	The mobile device is out of range.	Make sure the mobile device and the transmitter are located within 6 meters (20 feet). It is helpful to keep the devices on the same side of the body to minimize any radio frequency (RF) interference.
Lost communication alert	There is RF interference from other devices.	Move away from any equipment that can cause RF interference, such as cordless phones or routers.
Lost communication alert	The sensor disconnected from the transmitter.	Reconnect the sensor to the transmitter. Be careful not to pull the sensor out. Note: When the transmitter is reconnected to the sensor, the sensor will go through a warm-up period which may last up to 2 hours.

Problem	Likely Cause(s)	Resolution
Lost communication alert	The sensor pulled out from the skin.	Do not continue to use this sensor. Insert a new sensor to continue to receive SG values. For best results, recharge the transmitter before starting use of a new sensor to ensure full transmitter battery life. If the transmitter is still not communicating with the app, contact the local Medtronic representative for assistance.
Transmitter battery empty alert	The transmitter battery is empty and needs to be recharged.	See the <i>Guardian 4 Transmitter User Guide</i> for instructions on how to disconnect the transmitter from the sensor. For best results, recharge the transmitter between each use to ensure full transmitter battery life. Note: When the transmitter is reconnected to the sensor, the sensor will go through a warm-up period which may last up to 2 hours.
Mobile device battery low alert	The mobile device battery level is at 20% or lower. The battery needs to be recharged soon.	Recharge the mobile device battery to ensure that the Guardian 4 system can function and send alerts. Remember to always carry a charger for the mobile device to ensure continuous use of the Guardian 4 system.
Change sensor alert	The current sensor does not work properly and needs to be replaced.	To continue to receive SG values, a new sensor must be used. See the <i>Guardian 4 Transmitter User Guide</i> for instructions on how to change the sensor. For best results, recharge the transmitter between each use to ensure full transmitter battery life.
Sensor end of life alert	The current sensor has reached the end of its life and will no longer display SG values on the Guardian app.	To continue to receive SG values, a new sensor must be used. See the <i>Guardian 4 Transmitter User Guide</i> for instructions on how to change the sensor. For best results, recharge the transmitter between each use to ensure full transmitter battery life.
Calibration not accepted alert	The last calibration value entered was not accepted by the Guardian system.	Wait at least 15 minutes before attempting another calibration. The Guardian system will request another calibration, if needed, after 15 minutes from when the Calibration not accepted alert was received. Review the BG

Problem	Likely Cause(s)	Resolution
		meter instructions for use on how to test BG. Enter this new value in the app for calibration.
Calibrate now alert	A calibration is required by the system.	Review the BG meter instructions for use on how to test BG. Enter this new value in the app for calibration.
Sensor updating alert	There is an error with the sensor.	No action is required. The sensor is updating and this may take up to 2 hours. During this time, do not rely on alerts from the system as SG information will not be available. Monitor BG levels with the meter.
Transmitter error alert	There is an error with the transmitter.	Disconnect the transmitter from the sensor and reconnect it. Upon reconnecting the transmitter to the sensor, the sensor will go through a warm-up period which may last up to 2 hours. The transmitter may need to be replaced if issues are still being experienced. Contact 24-Hour Technical Support or the local Medtronic representative for more assistance.
Jailbroken device detected alert	The software on the iOS™ mobile device has been changed so that it no longer works in the way that the manufacturer intended.	The app cannot be used on a jailbroken device. To use the Guardian app the device software must remain as the manufacturer has designed. Changing the mobile device operating software causes the app to stop working.
Rooted device detected alert	The software on the Android™ mobile device has been changed so that it no longer works in the way that the manufacturer intended.	The app cannot be used on a rooted device. To use the Guardian app the device software must remain as the manufacturer has designed. Changing the mobile device operating software causes the app to stop working.

Maintenance

Cleaning

Not applicable.

Storage

Not applicable.

Disposal

Not applicable.

Technical specifications

Data security

The Guardian app has been designed with security features to help keep its data secure. However, there are important recommended steps to take to ensure the compatible mobile device used with the Guardian app is also secure:

- Do not leave the compatible mobile device unattended.
- Use caution when viewing or sharing data with others.
- Enable a security lock on the compatible mobile device. When the compatible mobile device is not in use, lock it in a way that requires the password to be entered in order to use it.
- Do not remove or interfere with the security features on the compatible mobile device.
- Do not attempt to modify the operating system, jailbreak, or root the device.
- Use only the official application store, such as the AppleTM* App StoreTM* or the Google PlayTM* store to get all mobile applications used with the compatible mobile device.
- Do not click on links from email messages, web pages, or text messages received from an unknown or untrusted source.
- Avoid the use of unknown Wi-FiTM* networks or public Wi-FiTM* hotspots.
- Enable security protection on your home Wi-FiTM* network, such as the use of a password and encryption.
- The app may send anonymous analytic data to Medtronic if permission has been granted in the setup of the app. This data is used to analyze crash logs and app performance. This access can be revoked or reinstated at any time in the Help screen of the app.

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USA:

Medtronic Diabetes Global
Headquarters
24-Hour Technical Support:
+1-800-646-4633
To order supplies:
+1-800-843-6687

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Introduction

The charger (MMT-7715) is an accessory intended for use with Medtronic continuous glucose monitoring (CGM) transmitter devices. Transmitters contain a non-replaceable, rechargeable battery. The charger is used to recharge a transmitter's internal battery.

Indications for use

The charger is indicated for use with Medtronic CGM transmitters. Refer to the system user guide for additional information, such as the intended age group.

Contraindications

No contraindications are associated with charger use. For contraindications related to continuous glucose monitoring, see the system user guide.

Clinical benefits

Not applicable. The charger is an accessory that is used to recharge the internal battery of Medtronic transmitters.

User safety

Warnings

Keep the charger out of reach of children. The charger contains small parts and may pose a choking hazard that can result in serious injury or death.

For questions or concerns related to charger use, contact a local Medtronic support representative for assistance.

For medical questions or concerns, or to report a serious incident, contact a healthcare professional.

Precautions

Consult the system user guide for precautions related to continuous glucose monitoring.

Risks and side effects

The charger contains small parts and may pose a choking hazard that can result in serious injury or death.

Hazardous substances

None.

Allergens

None known.

Using the charger

The charger has a green light that shows the charging status and a red light that communicates


any problems during charging. If you see a red light, see the **Troubleshooting** section. The charger requires one AAA alkaline battery.

Note: If the battery is installed incorrectly, or is low, the charger will not work. Repeat the battery installation steps using a new battery.

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To install a battery in the charger:

1. Push the battery cover in and slide it off, as shown in image A.
2. Insert a new AAA alkaline battery. Make sure the + and - symbols on the battery align with these same symbols shown on the charger.
3. Slide the cover back on the charger until it clicks into place, as shown in image B.

Troubleshooting

Problem	Likely Cause(s)	Resolution
A device is connected to the charger and no lights come on.	The charger battery has no power or no battery is inserted.	If there is no damage to the connector pins, replace the battery in the charger.
During charging, the flashing green light on the charger turns off and a longer flashing red light appears on the charger.	The charger battery is low on power.	Replace the battery in the charger.
During charging, a mix of quick and long flashing red lights appear on the charger.	The charger <i>and</i> the device are low on power.	Replace the battery in the charger and fully charge the device.

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To clean the charger:

1. Wash your hands thoroughly.
2. Use a damp cloth with mild cleaning solution, such as a dish washing detergent, to clean any dirt or foreign material from the outside of the charger. Never use organic solvents such as paint thinner or acetone to clean the charger.
3. Place the charger on a clean, dry cloth and air dry for two to three minutes.

Storage

Store the charger in a clean, dry location at room temperature.

Disposal


Dispose of the charger according to local regulations for electronic waste disposal.

Technical specifications

Operating Conditions	Temperature: 10° C to 40° C (50° F to 104° F) Relative Humidity: 30 % to 75 %, non-condensing
Storage Conditions	Temperature: -10° C to 50° C (14° F to 122° F) Relative Humidity: 10 % to 95 %, non-condensing

Icon Table

	Consult instructions for use
	Caution: consult instructions for use for important warnings or precautions not found on the label
	Batch code
	Catalogue number
	One per container/package

	Storage humidity limits
	Storage temperature limits
	Manufacturer
	Date of manufacture
	Medical device
	Importer
	Do not dispose of this product in unsorted municipal waste stream
	Authorized representative in the European Community
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts

© 2019 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic.

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ENGLISH

Welcome to the Guardian™ 4 system

Congratulations! Welcome to the Guardian™ 4 system community.

Access the continuous glucose monitoring system directly on a mobile device. With the Guardian™ 4 system, a mobile device can display sensor glucose information and provide alerts. To get started, follow the steps.

1

Charge the Guardian™ 4 transmitter

Place the transmitter on the charger. The green light on the charger will flash. The transmitter is fully charged when the green light turns off.

2

CareLink™

Create a CareLink™ account

Visit the CareLink™ website to create an account or log in with an existing CareLink™ Personal account.

3

Download the Guardian™ app

Search for and download the Guardian™ app from the Apple™ App Store™ or the Google Play™ store on the supported mobile device.

4

Learn how to use the Guardian™ 4 system

Use the information provided in the getting started guide, user guide, and Guardian™ app Help section to learn how to use the Guardian™ 4 system. Contact the local Medtronic support representative for more assistance.

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Medtronic

Guardian™ 4 System Performance

NOTE: You should review the information in this section with your healthcare professional.

NOTE: This document describes the performance of the Guardian 4 system (Guardian 4 sensor and Guardian 4 transmitter) when used with the Guardian App.

Clinical study description

The performance of the Guardian 4 system was evaluated using data collected during a multi-center prospective clinical study¹.

The study included participants ages 2 to 80 years old. The study enrolled a total of 335 subjects previously diagnosed with type 1 or 2 diabetes and 286 of these subjects completed the study. Subjects 18 and older were instructed to wear a total of three sensors and transmitters in the abdomen and arm. Subjects age 7 to 17 were instructed to wear a total of three sensors and transmitters in the arm and buttock. Subjects age 2 to 6 were instructed to wear a total of two sensors and transmitters in the arm and buttock. For all subjects, the transmitters were used to record raw sensor signals during the study and there was no real-time calculation of sensor glucose values.

Frequent sample testing (FST) was performed on four occasions for subjects 14 and older and on two occasions for subjects 2 to 13 years of age.

Reference blood (plasma) glucose values were obtained with a Yellow Springs Instrument (YSI®) Glucose Analyzer every 5-15 minutes for subjects 7 and older. Fingerstick blood glucose meter measurements were used as reference values for pediatric subjects aged 2 to 6 years old. During each FST, subjects with an established insulin sensitivity ratio and insulin carbohydrate ratio underwent a hypoglycemic challenge and a hyperglycemic challenge.

Data collected during the study was post-processed after the study using the Guardian 4 System sensor algorithm to convert the raw sensor information to sensor glucose values every five minutes. For the accuracy information presented in the following section, YSI reference values were paired with the closest sensor glucose reading within five minutes of the time of the reference value measurement.

NOTE: The Guardian 4 system is indicated for ages 7 and above. The data presented in this document is for ages 7 to 80 only.

Sensor accuracy

The overall sensor accuracy was calculated for sensors compared to a YSI® reference for 7 and older.

Table 1. Overall Accuracy Compared to YSI.

Patient Population	Insertion Site	Number of Subjects	Number of paired SG-YSI	Percent of SG within 20/20% of YSI	Mean Absolute Relative Difference (%)
Adults (18+)	Arm	153	20612	88.3 (88.0)	10.6
	Abdomen	150	18423	88.8 (88.4)	10.8
Pediatrics (7-17)	Arm	107	7702	85.6 (85.0)	11.6
	Buttocks	108	7624	82.7 (82.0)	12.3

* CGM readings are within 50-400 mg/dL, inclusive

*For 20% agreement, 20 mg/dL used when YSI<70 mg/dL

¹ Medtronic Inc., Performance Evaluation of an Advanced Algorithm with CGM in Adults, Adolescents and Pediatrics CIP324, September 2019

Table 2. Guardian 4 system accuracy, Adults, Arm insertion site.

Reference Range	Characteristic	FST 1	FST 2	FST 3	FST 4	Total
<3.9 mmol/L (<70 mg/dL)	No. of Paired Points	932	1021	535	440	2928
	Mean (SD)	17.60 (14.96)	10.43 (8.78)	11.06 (9.04)	8.97 (6.52)	12.61 (11.46)
	Median	13.9	8.45	8.02	7.88	9.47
3.9-10 mmol/L (70-180 mg/dL)	No. of Paired Points	3225	3037	2863	2363	11488
	Mean (SD)	12.85 (10.52)	9.63 (7.62)	9.69 (7.91)	7.86 (6.47)	10.19 (8.59)
	Median	10.54	8.05	7.69	6.42	8.06
>10 mmol/L (>180 mg/dL)	No. of Paired Points	1521	1446	1744	1485	6196
	Mean (SD)	13.71 (10.12)	10.51 (7.79)	9.21 (6.73)	8.86 (7.30)	10.53 (8.27)
	Median	11.94	9.09	8.04	7.54	8.75
Overall	No. of Paired Points	5678	5504	5142	4288	20612
	Mean (SD)	13.86 (11.39)	10.01 (7.90)	9.67 (7.68)	8.32 (6.79)	10.64 (9.00)
	Median	11.46	8.36	7.83	6.96	8.47

Table 3. Guardian 4 system accuracy, Adults, Abdomen insertion site.

Reference Range	Characteristic	FST 1	FST 2	FST 3	FST 4	Total
<3.9 mmol/L (<70 mg/dL)	No. of Paired Points	969	921	563	268	2721
	Mean (SD)	18.07 (15.35)	13.49 (10.65)	11.60 (10.56)	14.04 (9.48)	14.78 (12.67)
	Median	13.93	10.6	8.96	13.35	11.31
3.9-10 mmol/L (70-180 mg/dL)	No. of Paired Points	3244	2865	2504	1744	10357
	Mean (SD)	12.78 (10.42)	9.97 (7.93)	9.56 (7.91)	8.14 (7.18)	10.44 (8.84)
	Median	10.13	8.49	7.76	6.25	8.21
>10 mmol/L (>180 mg/dL)	No. of Paired Points	1418	1365	1366	1196	5345
	Mean (SD)	10.93 (8.14)	10.42 (7.31)	8.22 (6.55)	7.70 (6.36)	9.38 (7.29)
	Median	9.57	9.12	6.7	6.23	7.69
Overall	No. of Paired Points	5631	5151	4433	3208	18423
	Mean (SD)	13.23 (11.19)	10.72 (8.43)	9.41 (7.98)	8.47 (7.31)	10.78 (9.27)
	Median	10.49	8.92	7.49	6.5	8.39

Table 4. Guardian 4 system accuracy, Pediatrics, (7-17), Arm insertion site.

Reference Range	Characteristic	FST 1	FST 2	FST 3	FST 4	Total
<3.9 mmol/L (<70 mg/dL)	No. of Paired Points	293	255	103	81	732
	Mean (SD)	20.86 (16.67)	13.45 (9.89)	11.52 (6.98)	8.88 (6.29)	15.64 (13.27)
	Median	19.25	11.37	9.02	6.73	12.27

Reference Range	Characteristic	FST 1	FST 2	FST 3	FST 4	Total
3.9-10 mmol/L (70-180 mg/dL)	No. of Paired Points	1552	1167	386	384	3489
	Mean (SD)	11.86 (9.16)	12.13 (9.21)	11.49 (7.91)	10.96 (8.45)	11.81 (8.97)
	Median	10.26	10.26	10.21	9.16	10.12
>10 mmol/L (>180 mg/dL)	No. of Paired Points	1282	1124	656	419	3481
	Mean (SD)	12.01 (7.29)	10.36 (9.15)	10.47 (8.63)	7.49 (6.39)	10.64 (8.21)
	Median	11.45	8.17	8.01	5.53	9.19
Overall	No. of Paired Points	3127	2546	1145	884	7702
	Mean (SD)	12.77 (9.80)	11.48 (9.31)	10.91 (8.26)	9.12 (7.52)	11.65 (9.25)
	Median	11.27	9.43	8.81	7.21	9.77

Table 5. Guardian 4 system accuracy, Pediatrics, (7-17), Buttock insertion site.

Reference Range	Characteristic	FST 1	FST 2	FST 3	FST 4	Total
<3.9 mmol/L (<70 mg/dL)	No. of Paired Points	305	280	153	85	823
	Mean (SD)	21.33 (18.34)	11.79 (12.33)	19.12 (16.84)	16.38 (11.61)	17.16 (16.10)
	Median	15.47	7.48	15.19	13.49	11.71
3.9-10 mmol/L (70-180 mg/dL)	No. of Paired Points	1566	1167	451	309	3493
	Mean (SD)	12.52 (10.71)	10.89 (9.53)	13.39 (10.80)	13.48 (11.40)	12.17 (10.45)
	Median	9.94	8.28	11.94	10.24	9.58
>10 mmol/L (>180 mg/dL)	No. of Paired Points	1222	1061	670	355	3308
	Mean (SD)	12.14 (8.94)	10.32 (8.17)	11.28 (10.50)	9.83 (7.47)	11.13 (8.94)
	Median	10.16	8.29	8.22	8.7	8.97
Overall	No. of Paired Points	3093	2508	1274	749	7624
	Mean (SD)	13.24 (11.40)	10.75 (9.36)	12.97 (11.79)	12.08 (10.02)	12.26 (10.76)
	Median	10.4	8.18	9.88	9.77	9.5

Alert performance

Glucose True Alert Rate

The glucose true alert rate is the rate at which the blood glucose confirmed that the CGM alert was triggered correctly.

Table 6. Glucose True Alert Rate, Adults.

Glucose True Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Abdomen	31.6%	28.9%	18.3%	15.0%	20.3%	17.1%
	Arm	27.3%	25.8%	17.5%	14.6%	19.5%	16.8%

Glucose True Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
3 mmol/L (54 mg/dL)	Abdomen	50.0%	48.3%	28.0%	23.6%	32.1%	28.2%
	Arm	51.6%	51.6%	26.0%	22.3%	31.8%	28.9%
3.3 mmol/L (60 mg/dL)	Abdomen	75.9%	75.0%	43.4%	40.9%	51.6%	49.5%
	Arm	73.3%	71.9%	42.0%	38.9%	50.7%	48.0%
3.9 mmol/L (70 mg/dL)	Abdomen	86.9%	86.3%	54.6%	51.6%	64.5%	62.2%
	Arm	85.9%	85.9%	51.4%	48.4%	62.2%	60.2%
4.4 mmol/L (80 mg/dL)	Abdomen	85.0%	84.1%	54.5%	51.8%	64.5%	62.4%
	Arm	81.9%	80.8%	49.0%	46.4%	59.7%	57.6%
5 mmol/L (90 mg/dL)	Abdomen	80.0%	80.0%	51.5%	48.7%	61.2%	59.4%
	Arm	77.7%	77.4%	46.1%	43.1%	56.3%	54.2%
10 mmol/L (180 mg/dL)	Abdomen	85.5%	85.3%	61.7%	59.4%	70.5%	68.9%
	Arm	87.6%	87.2%	66.4%	64.0%	74.3%	72.7%
12.2 mmol/L (220 mg/dL)	Abdomen	83.8%	83.2%	60.6%	58.0%	68.9%	67.1%
	Arm	86.8%	86.8%	64.2%	62.1%	72.2%	70.9%
13.9 mmol/L (250 mg/dL)	Abdomen	85.5%	85.0%	60.9%	58.1%	69.3%	67.3%
	Arm	87.9%	87.0%	67.3%	62.0%	74.5%	70.8%
16.7 mmol/L (300 mg/dL)	Abdomen	87.1%	86.0%	64.9%	59.8%	72.7%	68.9%
	Arm	88.9%	88.9%	69.5%	65.5%	76.1%	73.5%

Table 7. Glucose True Alert Rate, Pediatrics, (7-17).

Glucose True Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Buttock	20.0%	20.0%	10.1%	8.9%	12.0%	11.0%
	Arm	23.8%	23.8%	9.8%	7.7%	13.0%	11.4%
3 mmol/L (54 mg/dL)	Buttock	29.6%	29.6%	13.1%	11.9%	17.0%	16.1%
	Arm	32.0%	32.0%	12.1%	10.9%	16.5%	15.6%
3.3 mmol/L (60 mg/dL)	Buttock	60.5%	60.5%	30.0%	26.8%	38.1%	35.6%
	Arm	52.8%	52.8%	30.4%	27.2%	36.5%	34.2%
3.9 mmol/L (70 mg/dL)	Buttock	71.9%	70.8%	36.8%	33.2%	46.6%	43.7%
	Arm	63.0%	63.0%	35.7%	32.2%	43.4%	40.8%
4.4 mmol/L (80 mg/dL)	Buttock	68.1%	67.4%	40.8%	39.1%	49.9%	48.5%
	Arm	69.1%	69.1%	42.4%	39.2%	51.2%	49.1%
5 mmol/L (90 mg/dL)	Buttock	71.0%	71.0%	49.0%	47.1%	57.2%	56.0%
	Arm	74.4%	72.2%	48.4%	44.2%	57.8%	54.3%

Glucose True Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
10 mmol/L (180 mg/dL)	Buttock	91.4%	90.9%	76.4%	74.0%	82.7%	81.0%
	Arm	92.2%	92.2%	82.7%	78.3%	86.8%	84.3%
12.2 mmol/L (220 mg/dL)	Buttock	90.5%	90.5%	74.3%	72.4%	80.7%	79.6%
	Arm	91.3%	90.2%	77.3%	74.0%	83.0%	80.5%
13.9 mmol/L (250 mg/dL)	Buttock	88.4%	87.8%	71.5%	68.0%	77.9%	75.5%
	Arm	88.4%	88.4%	70.7%	68.4%	77.4%	75.9%
16.7 mmol/L (300 mg/dL)	Buttock	79.2%	79.2%	54.4%	48.7%	62.6%	58.7%
	Arm	85.2%	85.2%	60.5%	59.9%	68.9%	68.5%

Glucose False Alert Rate

The glucose false alert rate is the rate at which the BG did not confirm that the CGM alert was triggered correctly.

Table 8. Glucose False Alert Rate, Adults.

Glucose False Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Abdomen	68.4%	71.1%	81.7%	85.0%	79.7%	82.9%
	Arm	72.7%	74.2%	82.5%	85.4%	80.5%	83.2%
3 mmol/L (54 mg/dL)	Abdomen	50.0%	51.7%	72.0%	76.4%	67.9%	71.8%
	Arm	48.4%	48.4%	74.0%	77.7%	68.2%	71.1%
3.3 mmol/L (60 mg/dL)	Abdomen	24.1%	25.0%	56.6%	59.1%	48.4%	50.5%
	Arm	26.7%	28.1%	58.0%	61.1%	49.3%	52.0%
3.9 mmol/L (70 mg/dL)	Abdomen	13.1%	13.7%	45.4%	48.4%	35.5%	37.8%
	Arm	14.1%	14.1%	48.6%	51.6%	37.8%	39.8%
4.4 mmol/L (80 mg/dL)	Abdomen	15.0%	15.9%	45.5%	48.2%	35.5%	37.6%
	Arm	18.1%	19.2%	51.0%	53.6%	40.3%	42.4%
5 mmol/L (90 mg/dL)	Abdomen	20.0%	20.0%	48.5%	51.3%	38.8%	40.6%
	Arm	22.3%	22.6%	53.9%	56.9%	43.7%	45.8%
10 mmol/L (180 mg/dL)	Abdomen	14.5%	14.7%	38.3%	40.6%	29.5%	31.1%
	Arm	12.4%	12.8%	33.6%	36.0%	25.7%	27.3%
12.2 mmol/L (220 mg/dL)	Abdomen	16.2%	16.8%	39.4%	42.0%	31.1%	32.9%
	Arm	13.2%	13.2%	35.8%	37.9%	27.8%	29.1%
13.9 mmol/L (250 mg/dL)	Abdomen	14.5%	15.0%	39.1%	41.9%	30.7%	32.7%
	Arm	12.1%	13.0%	32.7%	38.0%	25.5%	29.2%

Glucose False Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
16.7 mmol/L (300 mg/dL)	Abdomen	12.9%	14.0%	35.1%	40.2%	27.3%	31.1%
	Arm	11.1%	11.1%	30.5%	34.5%	23.9%	26.5%

Table 9. Glucose False Alert Rate, Pediatrics, (7-17).

Glucose False Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Buttock	80.0%	80.0%	89.9%	91.1%	88.0%	89.0%
	Arm	76.2%	76.2%	90.2%	92.3%	87.0%	88.6%
3 mmol/L (54 mg/dL)	Buttock	70.4%	70.4%	86.9%	88.1%	83.0%	83.9%
	Arm	68.0%	68.0%	87.9%	89.1%	83.5%	84.4%
3.3 mmol/L (60 mg/dL)	Buttock	39.5%	39.5%	70.0%	73.2%	61.9%	64.4%
	Arm	47.2%	47.2%	69.6%	72.8%	63.5%	65.8%
3.9 mmol/L (70 mg/dL)	Buttock	28.1%	29.2%	63.2%	66.8%	53.4%	56.3%
	Arm	37.0%	37.0%	64.3%	67.8%	56.6%	59.2%
4.4 mmol/L (80 mg/dL)	Buttock	31.9%	32.6%	59.2%	60.9%	50.1%	51.5%
	Arm	30.9%	30.9%	57.6%	60.8%	48.8%	50.9%
5 mmol/L (90 mg/dL)	Buttock	29.0%	29.0%	51.0%	52.9%	42.8%	44.0%
	Arm	25.6%	27.8%	51.6%	55.8%	42.2%	45.7%
10 mmol/L (180 mg/dL)	Buttock	8.6%	9.1%	23.6%	26.0%	17.3%	19.0%
	Arm	7.8%	7.8%	17.3%	21.7%	13.2%	15.7%
12.2 mmol/L (220 mg/dL)	Buttock	9.5%	9.5%	25.7%	27.6%	19.3%	20.4%
	Arm	8.7%	9.8%	22.7%	26.0%	17.0%	19.5%
13.9 mmol/L (250 mg/dL)	Buttock	11.6%	12.2%	28.5%	32.0%	22.1%	24.5%
	Arm	11.6%	11.6%	29.3%	31.6%	22.6%	24.1%
16.7 mmol/L (300 mg/dL)	Buttock	20.8%	20.8%	45.6%	51.3%	37.4%	41.3%
	Arm	14.8%	14.8%	39.5%	40.1%	31.1%	31.5%

Glucose Correct Detection Alert Rate

Glucose Correct Detection Rate is the rate that the device alerted when it should have alerted.

Table 10. Glucose Correct Detection Alert Rate, Adults.

Glucose Correct Detection Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Abdomen	27.3%	25.0%	81.8%	72.7%	81.8%	72.7%
	Arm	37.5%	35.4%	91.7%	83.3%	91.7%	85.4%
3 mmol/L (54 mg/dL)	Abdomen	40.0%	37.3%	85.3%	77.3%	86.7%	77.3%
	Arm	56.3%	56.3%	87.4%	82.8%	87.4%	83.9%
3.3 mmol/L (60 mg/dL)	Abdomen	55.1%	54.5%	89.7%	85.9%	89.7%	86.5%
	Arm	65.1%	64.0%	91.3%	88.4%	91.9%	89.0%
3.9 mmol/L (70 mg/dL)	Abdomen	78.0%	78.0%	97.8%	95.2%	97.8%	96.2%
	Arm	87.3%	86.8%	97.2%	96.2%	97.6%	96.7%
4.4 mmol/L (80 mg/dL)	Abdomen	83.3%	82.1%	94.4%	91.5%	94.4%	92.3%
	Arm	84.2%	83.1%	92.5%	89.5%	92.5%	90.2%
5 mmol/L (90 mg/dL)	Abdomen	82.1%	81.8%	93.8%	89.7%	93.8%	90.7%
	Arm	87.3%	86.3%	96.2%	92.1%	96.2%	93.0%
10 mmol/L (180 mg/dL)	Abdomen	87.2%	86.3%	95.9%	92.3%	95.9%	94.1%
	Arm	83.0%	81.8%	94.2%	90.4%	94.4%	91.3%
12.2 mmol/L (220 mg/dL)	Abdomen	85.4%	85.1%	95.7%	93.6%	95.7%	94.2%
	Arm	79.4%	78.3%	92.4%	90.6%	92.7%	90.9%
13.9 mmol/L (250 mg/dL)	Abdomen	76.1%	75.6%	90.2%	87.6%	90.2%	88.0%
	Arm	72.5%	71.1%	88.3%	85.3%	88.3%	86.1%
16.7 mmol/L (300 mg/dL)	Abdomen	73.5%	71.7%	90.3%	85.0%	91.2%	85.8%
	Arm	62.3%	60.1%	83.3%	80.4%	83.3%	80.4%

Table 11. Glucose Correct Detection Alert Rate, Pediatrics, (7-17).

Glucose Correct Detection Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Buttock	50.0%	50.0%	93.8%	87.5%	93.8%	87.5%
	Arm	66.7%	66.7%	86.7%	73.3%	86.7%	80.0%
3 mmol/L (54 mg/dL)	Buttock	69.6%	69.6%	91.3%	87.0%	91.3%	87.0%
	Arm	72.7%	72.7%	90.9%	86.4%	90.9%	86.4%
3.3 mmol/L (60 mg/dL)	Buttock	71.0%	69.6%	95.7%	85.5%	95.7%	88.4%
	Arm	64.6%	64.6%	89.2%	83.1%	89.2%	84.6%
3.9 mmol/L (70 mg/dL)	Buttock	83.5%	81.2%	98.8%	92.9%	98.8%	96.5%
	Arm	83.3%	81.0%	96.4%	90.5%	96.4%	95.2%

Glucose Correct Detection Alert Rate

Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
4.4 mmol/L (80 mg/dL)	Buttock	87.0%	85.2%	96.5%	93.9%	97.4%	95.7%
	Arm	86.6%	85.7%	97.3%	93.8%	97.3%	95.5%
5 mmol/L (90 mg/dL)	Buttock	86.1%	84.2%	98.7%	94.9%	99.4%	95.6%
	Arm	89.9%	87.9%	98.0%	91.9%	98.7%	94.6%
10 mmol/L (180 mg/dL)	Buttock	91.5%	90.9%	97.1%	92.0%	97.1%	95.5%
	Arm	90.3%	89.0%	96.2%	91.3%	96.9%	94.1%
12.2 mmol/L (220 mg/dL)	Buttock	86.0%	83.1%	93.4%	90.4%	93.8%	91.9%
	Arm	88.0%	84.9%	95.9%	92.1%	96.2%	94.5%
13.9 mmol/L (250 mg/dL)	Buttock	78.1%	77.6%	92.9%	88.3%	93.4%	89.8%
	Arm	77.1%	74.8%	91.6%	85.0%	91.6%	87.4%
16.7 mmol/L (300 mg/dL)	Buttock	65.3%	65.3%	84.2%	76.8%	85.3%	78.9%
	Arm	72.2%	72.2%	93.8%	88.7%	93.8%	89.7%

Glucose Missed Detection Alert Rate

Glucose Missed Detection Rate is the rate that the device did not alert when it should have alerted.

Table 12. Glucose Missed Detection Alert Rate, Adults.

Glucose Missed Detection Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Abdomen	72.7%	75.0%	18.2%	27.3%	18.2%	27.3%
	Arm	62.5%	64.6%	8.3%	16.7%	8.3%	14.6%
3 mmol/L (54 mg/dL)	Abdomen	60.0%	62.7%	14.7%	22.7%	13.3%	22.7%
	Arm	43.7%	43.7%	12.6%	17.2%	12.6%	16.1%
3.3 mmol/L (60 mg/dL)	Abdomen	44.9%	45.5%	10.3%	14.1%	10.3%	13.5%
	Arm	34.9%	36.0%	8.7%	11.6%	8.1%	11.0%
3.9 mmol/L (70 mg/dL)	Abdomen	22.0%	22.0%	2.2%	4.8%	2.2%	3.8%
	Arm	12.7%	13.2%	2.8%	3.8%	2.4%	3.3%
4.4 mmol/L (80 mg/dL)	Abdomen	16.7%	17.9%	5.6%	8.5%	5.6%	7.7%
	Arm	15.8%	16.9%	7.5%	10.5%	7.5%	9.8%
5 mmol/L (90 mg/dL)	Abdomen	17.9%	18.2%	6.2%	10.3%	6.2%	9.3%
	Arm	12.7%	13.7%	3.8%	7.9%	3.8%	7.0%
10 mmol/L (180 mg/dL)	Abdomen	12.8%	13.7%	4.1%	7.7%	4.1%	5.9%
	Arm	17.0%	18.2%	5.8%	9.6%	5.6%	8.7%
12.2 mmol/L (220 mg/dL)	Abdomen	14.6%	14.9%	4.3%	6.4%	4.3%	5.8%
	Arm	20.6%	21.7%	7.6%	9.4%	7.3%	9.1%

Glucose Missed Detection Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
13.9 mmol/L (250 mg/dL)	Abdomen	23.9%	24.4%	9.8%	12.4%	9.8%	12.0%
	Arm	27.5%	28.9%	11.7%	14.7%	11.7%	13.9%
16.7 mmol/L (300 mg/dL)	Abdomen	26.5%	28.3%	9.7%	15.0%	8.8%	14.2%
	Arm	37.7%	39.9%	16.7%	19.6%	16.7%	19.6%

Table 13. Glucose Missed Detection Alert Rate, Pediatrics, (7-17).

Glucose Missed Detection Alert Rate							
Glucose	Insertion Site	Threshold Only		Predictive Only		Threshold & Predictive	
		30 min	15 min	30 min	15 min	30 min	15 min
2.8 mmol/L (50 mg/dL)	Buttocks	50.0%	50.0%	6.3%	12.5%	6.3%	12.5%
	Arm	33.3%	33.3%	13.3%	26.7%	13.3%	20.0%
3 mmol/L (54 mg/dL)	Buttocks	30.4%	30.4%	8.7%	13.0%	8.7%	13.0%
	Arm	27.3%	27.3%	9.1%	13.6%	9.1%	13.6%
3.3 mmol/L (60 mg/dL)	Buttock	29.0%	30.4%	4.3%	14.5%	4.3%	11.6%
	Arm	35.4%	35.4%	10.8%	16.9%	10.8%	15.4%
3.9 mmol/L (70 mg/dL)	Buttock	16.5%	18.8%	1.2%	7.1%	1.2%	3.5%
	Arm	16.7%	19.0%	3.6%	9.5%	3.6%	4.8%
4.4 mmol/L (80 mg/dL)	Buttock	13.0%	14.8%	3.5%	6.1%	2.6%	4.3%
	Arm	13.4%	14.3%	2.7%	6.3%	2.7%	4.5%
5 mmol/L (90 mg/dL)	Buttock	13.9%	15.8%	1.3%	5.1%	0.6%	4.4%
	Arm	10.1%	12.1%	2.0%	8.1%	1.3%	5.4%
10 mmol/L (180 mg/dL)	Buttock	8.5%	9.1%	2.9%	8.0%	2.9%	4.5%
	Arm	9.7%	11.0%	3.8%	8.7%	3.1%	5.9%
12.2 mmol/L (220 mg/dL)	Buttock	14.0%	16.9%	6.6%	9.6%	6.3%	8.1%
	Arm	12.0%	15.1%	4.1%	7.9%	3.8%	5.5%
13.9 mmol/L (250 mg/dL)	Buttock	21.9%	22.4%	7.1%	11.7%	6.6%	10.2%
	Arm	22.9%	25.2%	8.4%	15.0%	8.4%	12.6%
16.7 mmol/L (300 mg/dL)	Buttock	34.7%	34.7%	15.8%	23.2%	14.7%	21.1%
	Arm	27.8%	27.8%	6.2%	11.3%	6.2%	10.3%

Sensor life

Adults

95.7% of sensors worn in the abdomen operated more than six days and up to the full seven days of wear (144 to 168 hours). The mean functional sensor life for sensors worn in the abdomen insertion site over the course of the study was 138.6 hours, with a median functional life of 164.7 hours. 93.4% of sensors

worn in the arm operated more than six days and up to the full seven days of wear (144 to 168 hours). The mean functional sensor life for sensors worn in the arm insertion site over the course of the study was 147.9 hours, with a median functional life of 167.9 hours.

Pediatrics 7-17

95.4% of sensors worn in the arm operated more than six days and up to the full seven days of wear (144 to 168 hours). The mean functional sensor life for sensors worn in the arm insertion site over the course of the study was 138.9 hours, with a median functional life of 163.7 hours. 93.5% of sensors worn in the buttock operated more than six days and up to the full seven days of wear (144 to 168 hours). The mean functional sensor life for sensors worn in the buttock insertion site over the course of the study was 137.2 hours, with a median functional life of 152.8 hours.

Safety

Device related adverse events were limited to mild bleeding or irritation at the sensor insertion site or syncopal or vasovagal episodes related to sensor insertion.

© 2021 Medtronic. All rights reserved.

Medtronic, Medtronic logo and Further, Together are trademarks of Medtronic. Third party brands are trademarks of their respective owners. All other brands are trademarks of a Medtronic company

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United States

FCC ID: OH27821

This device complies with part 15 of the FCC Rules. Operation is subject to two conditions:

(1) This device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This device has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation.

This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this device does cause harmful interference to radio or television reception, which can be determined by turning the device off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna.
- Increase the separation between the device and the receiver.
- Decrease the distance between the transmitter and the compatible mobile device to 6.1 meters (20 feet) or less.
- Increase the separation between the transmitter and the equipment that is receiving/emitting interference.

IMPORTANT: Do not change or modify the internal RF transmitter or antenna unless expressly approved by Medtronic Diabetes. Doing so could interfere with your ability to operate the device.

Australia and New Zealand

This product conforms to Australia and New Zealand Radio Requirements.

Canada

IC:3408B-7821

This device contains license-exempt transmitter(s)/receiver(s) that comply with Innovation, Science and Economic Development Canada's license-exempt RSS(s). Operation is subject to the following two conditions:

(1) this device may not cause interference, and (2) this device must accept any interference, including interference that may cause undesired operation of the device.

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Mexico

The operation of this equipment is subject to the following two conditions: (1) it is possible that this equipment or device may not cause harmful interference and (2) this equipment or device must accept any interference, including interference that may cause its undesired operation.

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Singapore

Complies with
IMDA Standards
DA103452

This product complies with IMDA Standards.
该产品符合 IMDA 标准

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07642-19-02516

This equipment is not entitled to protection against harmful interference and may not cause interference with duly authorized systems.

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(European Conformity)

EN Hereby, Medtronic declares that this Radio Equipment is in compliance with Directive 2014/53/EU.

The full text of the EU declaration of conformity is available at the following internet address:

www.MedtronicCompliance.eu

This product operates in the 2.4–2.4835 GHz band with a radiated output power of less than 0.97 mW e.i.r.p.

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24-Stunden-Hotline: 0800 633333

Serbia:

Epsilon Research International d.o.o.

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24 – Stunden – Hotline: 0820 820 190

Introduction

The tester (MMT-7736L) is an accessory item included with select Medtronic Diabetes continuous glucose transmitters.

Indications for use

The tester is used to test the transmitter. The tester is also used to form a waterproof seal that protects the transmitter connector while the transmitter is being cleaned.

Contraindications

None known.

Clinical benefits

Not applicable. The tester is an accessory that tests the transmitter and protects it during cleaning.

User safety

Warnings

Keep the tester out of reach of children. The tester may pose a choking hazard that can result in serious injury or death.

Do not use the tester if it comes into contact with blood. Contact with blood may cause infection. If the tester comes into contact with blood, dispose directly into a sharps container.

Precautions

Read the entire transmitter user guide to prevent misuse of the tester. Using the tester before reading the transmitter user guide can lead to misuse of the transmitter.

Risks and side effects

The tester may pose a choking hazard to children, which can result in serious injury or death.

Hazardous substances

None.

Allergens

Internal components contain nickel.

Using the tester

Refer to the transmitter user guide for instructions.

Maintenance

Cleaning

Not applicable.

Storage

Store the tester at room temperature.

Disposal


Dispose of the tester according to local regulations for electronic waste disposal.

Technical specifications

The tester has an expected lifetime of one year when connected to the transmitter once per sensor replacement period. If physical damage or discoloration is observed, the tester should be replaced.

Approximate dimensions	Approximate weight
3.8 x 2.0 x 1.0 centimeters (1.5 x 0.8 x 0.3 inches)	3.0 grams (0.11 ounces)

Icon table

Symbol	Meaning
	Manufacturer
	Authorized representative in the European Community.
	Catalogue number
	Storage temperature limits
	Consult instructions for use
	Medical device
	Importer
	Date of manufacture
(2x)	Two per container/package
	Batch code
	Storage humidity limits
	Do not dispose of this product in unsorted municipal waste stream
	Caution
	Conformité Européenne (European Conformity). This symbol means that the device fully complies with applicable European Union Acts.

© 2019 Medtronic. All rights reserved. Medtronic, Medtronic logo, and Further, Together are trademarks of Medtronic. All other brands are trademarks of a Medtronic company.

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USA:

Medtronic Diabetes Global
Headquarters
24-Hour Technical Support:
+1-800-646-4633
To order supplies: +1-800-843-6687

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
Introduction

The one-pressserter (MMT-7512) is used as an aid to insert Medtronic MiniMed™ glucose sensors.

- A. Bump (on both buttons)
- B. Thumbprint marking

Indications for use

The sserter is indicated for single-patient use. It is not intended for use by multiple patients. Refer to the system user guide for additional information, such as the intended age group.

Contraindications

No contraindications are associated with sserter use. For contraindications related to continuous glucose monitoring, see the system user guide.

Clinical benefits

Not applicable. The sserter is an accessory that helps to insert compatible Medtronic sensors.

User Safety

Warnings

Read this entire user guide prior to sensor insertion. The sserter does not work the same way as other Medtronic insertion devices. Failure to follow directions may result in improper insertion, pain, or injury.

For additional warnings related to sensor use, see the sensor user guide.

Keep the sserter out of reach of children. The sserter contains small parts and may pose a choking hazard that can result in serious injury or death.

Do not make therapy decisions based on sensor glucose values. Sensor glucose and blood glucose values may differ. If a sensor glucose reading is low or high, or symptoms of low or high blood glucose appear, use a BG meter to confirm blood glucose prior to making therapy decisions. Failure to check blood glucose levels with a BG meter prior to making therapy decisions can result in the infusion of too much or too little insulin, which can cause hypoglycemia or hyperglycemia.

Healthcare professionals and caregivers:

- A retractable needle is attached to the sensor. Minimal bleeding may occur. Always wear gloves to insert the sensor
- Cover the sensor with sterile gauze to remove the needle housing from the sensor

Place the needle housing directly into a sharps container after sensor insertion to prevent accidental needle stick injury.

Always inspect the sensor packaging for damage before use. Sensors are sterile and non-pyrogenic unless the package has been opened or damaged. If the sensor packaging is open or damaged, discard the sensor directly into a medical waste container. Use of a non-sterile sensor can result in infection at the insertion site.

For questions or concerns related to serter use, contact a local Medtronic support representative for assistance.
For medical questions or concerns, contact a healthcare provider.

Precautions

Wash hands with soap and water prior to sensor insertion to reduce the risk of infection.

Inspect the insertion site for bleeding under, around, or on top of the sensor. If bleeding occurs, use a sterile gauze pad to apply steady pressure to the sensor for three minutes.

WARNING: Do not use non-sterile gauze to apply pressure to the sensor. Use of non-sterile gauze can cause infection at the insertion site.

If bleeding stops after three minutes of steady pressure, connect the transmitter to the sensor.

CAUTION: Do not connect the transmitter to the sensor if bleeding at the insertion site does not stop. Active bleeding can damage the transmitter and result in a loss of communication between the transmitter and the receiving device.

If bleeding does not stop after three minutes of steady pressure, remove the sensor, and continue to apply pressure until bleeding stops. Discard the used sensor directly in a sharps container after removal. Insert a new sensor at a different insertion site.

Note: Always use a new insertion site to insert a new sensor. Reuse of the same insertion site can cause redness, irritation, pain, bleeding, tenderness, and inflammation.

Do not clean, resterilize, or attempt to remove the needle from the needle housing. Needle stick injury may occur.

Do not reuse the sensor. Reuse of the sensor can cause damage to the sensor surface and lead to inaccurate sensor glucose values, irritation, and infection at the insertion site.

Rotate the sensor insertion site so that sites do not become overused.

Confirm that the sensor is firmly attached to the skin to prevent it from coming out of the insertion site. Factors that increase the likelihood of the sensor coming out include increased physical activity (particularly in younger patients) and improper taping technique.

Risks and side effects

The serter contains small parts and may pose a choking hazard that can result in serious injury or death.

Side effects include discomfort and skin irritation at the insertion site.

Hazardous substances

None.

Allergens

None known.

Where to insert the sensor

Choose an insertion site with a sufficient layer of subcutaneous fat.

Note: Do not insert the sensor within 5.0 cm (2 inches) of the navel area to avoid discomfort at the insertion site and incomplete sensor adhesion.

For optimal sensor performance, and to prevent accidental sensor removal:

- Do not insert the sensor into muscle, tough skin, or scar tissue
- Avoid areas that are constrained by clothing or accessories
- Avoid areas that are located under a belt or waistband

- Avoid areas subjected to vigorous movement during exercise

For a diagram of recommended sensor insertion areas, see the sensor user guide.

Inserting the sensor

For instructions on how to insert the sensor, see the sensor user guide.

Maintenance

Cleaning

Clean the serter after every use or whenever debris, blood, or other contaminants are visible.

To clean the serter:

Note: Press and release the mechanism on the underside of the serter to fully rinse the device.

1. Rinse the serter under running tap water.
2. Submerge the serter in a solution of mild soap and warm water.
3. Use a soft-bristled brush to clean the serter until it is free of visible dirt.
4. Rinse the serter under running tap water until all visible soap is gone. Make sure to fully rinse all hard-to-reach areas.
5. Shake off any excess water. Place the serter upright on a clean cloth to air dry.

Storage

Store the serter in the released position at room temperature.



Disposal

Dispose of the serter according to local regulations.

Technical specifications

Approximate Dimensions	Approximate Weight
7.85 x 6.91 x 5.77 centimeters (3.09 x 2.72 x 2.27 inches)	112.5 grams (3.97 ounces)

Icon table

	Consult instructions for use
	Caution
	Batch code
	Catalogue number
	One per container/package
	Storage humidity limit
	Storage temperature limit
	Manufacturer
	Date of manufacture
	Importer
	Signifies European technical conformity
	Medical device
	Authorized representative in the European Community

©2019, Medtronic MiniMed, Inc. All rights reserved. MiniMed™ is a trademark of Medtronic MiniMed, Inc.

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State of California)

County of Ventura)


On March 7, 2023 before me, Shacoya Braggs, Notary Public

Date Here Insert Name and Title of the Officer
personally appeared Dolphine Ameriya Katiwa, PhD
Name(s) of Signer(s)

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Description of Attached Document

Title or Type of Document: Attestation Statement Document Date: March 7, 2023

Number of Pages: _____ Signer(s) Other Than Named Above: _____

Capacity(ies) Claimed by Signer(s)

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Signer's Name: _____

☐ Corporate Officer — Title(s): _____

☐ Partner — ☐ Limited ☐ General

☐ Individual ☐ Attorney in Fact

☐ Trustee ☐ Guardian or Conservator

☐ Other: _____

Signer Is Representing: _____

Medtronic

18000 Devonshire Street
Northridge, CA 91325 USA
www.medtronic.com
Tel 818-576-5211

ATTESTATION STATEMENT

As a responsible member of Medtronic MiniMed, I hereby attest that I believe the following Appendix To The Instruction For Use, attached herein, to be a true and accurate copy.

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- EN © 2021 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks ("™") belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.
© 2019 Medtronic. All rights reserved. Medtronic, Medtronic 로고 및 Further, Together는 Medtronic의 상표입니다. 기타 모든 브랜드는 Medtronic 기업의 상표입니다.