convatec
— forever caring —

Unomedical A/S
(a Convatec company)
Aaholmvej 1-3, Østød
4320 Lejre, Denmark
convatec.com
infusion-set.com

4 April 2024

Attestation statement

As a responsible member of Unomedical A/S (a Convatec company representing the Convatec Infusion Care business unit) with offices at Aaholmvej 1-3, Østød, 4320 Lejre, Denmark, I hereby attest, that I believe the following document: "Addendum to the instructions for use" attached herein, to be true and accurate copy.

Respectfully,
Lukas Solanko
Regulatory Affairs Specialist

Nina Schwarz
Regulatory Affairs Coordinator

Unomedical a/s
A ConvaTec Company
Aaholmvej 1-3, Østød
DK-4320 Lejre, Denmark
Reg./CVR-No: 6415 3315

Unomedical A/S (a Convatec company) CVR 64153315
Aaholmvej 1-3, Østød, 4320 Lejre, Denmark

© Convatec 2023

This is to certify that **Mrs Nina Schwarz** today in my presence at the notary office has approved and signed this document. No conspicuous corrections or addenda were found in the document.

A Danish notary has not reviewed or approved the contents of the document and have only verified the identity, the accuracy of the dating and the authenticity of the signature.

Mrs Nina Schwarz has proved her identity by presenting Danish Driver's License No 33622337.

District Court of Lyngby, Denmark, on Thursday, April 4, 2024

Ritha Ahle
Ritha Ahle
Notary Public

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- Insertion angle: 90 degrees
- 29 gauge needle
- 6mm, 8mm and 10mm stainless steel needle
- Tubing length: 46 cm (18"), 60 cm (23") and 80 cm (32")
- Tubing ID: 0.385 mm
- OD: 1.50 mm
- Usage time: Max 48 hours

MMT-862	SureT G29 6 mm, 18", Paradigm
MMT-864	SureT G29 6 mm, 23", Paradigm
MMT-866	SureT G29 6 mm, 32", Paradigm
MMT-874	SureT G29 8 mm, 23", Paradigm
MMT-876	SureT G29 8 mm, 32", Paradigm
MMT-884	SureT G29 10 mm, 23", Paradigm
MMT-886	SureT G29 10 mm, 32", Paradigm
MMT-864T	SureT G29 6 mm, 23", Paradigm
MMT-862A	SureT G29 6 mm, 18", (SC1)
MMT-864A	SureT G29 6 mm, 23", (SC1)
MMT-866A	SureT G29 6 mm, 32", (SC1)
MMT-874A	SureT G29 8 mm, 23", (SC1)
MMT-876A	SureT G29 8 mm, 32", (SC1)

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- Insertion angle: 90 degrees
- 27 gauge introducer needle
- Soft cannula OD 0.68 mm
- Soft cannula length: 6 mm and 9 mm
- Tubing length: 46 cm (18") 60 cm (23") and 80 cm (32")
- Tubing ID: 0.385 mm, OD: 1.50 mm
- Usage time: Max 72 hours

MMT-921	Pink, Paradigm, 6mm, 18"/46cm
MMT-923	Pink, Paradigm, 6mm, 23"/60cm
MMT-925	Pink, Paradigm, 6mm, 32"/80cm
MMT-941	Blue, Paradigm, 6mm, 18"/46cm
MMT-943	Blue, Paradigm, 6mm, 23"/60cm
MMT-945	Blue, Paradigm, 6mm, 32"/80cm
MMT-961	Clear, Paradigm, 6mm, 18"/46cm
MMT-963	Clear, Paradigm, 6mm, 23"/60cm
MMT-965	Clear, Paradigm, 6mm, 32"/80cm
MMT-975	Clear, Paradigm, 9mm, 32"/80cm
MMT-921A	Pink, SC1 6mm, 18"/46cm
MMT-923A	Pink, SC1, 6mm, 23"/60cm
MMT-925A	Pink, SC1 6mm, 32"/80cm
MMT-941A	Blue, SC1, 6mm, 18"/46cm
MMT-943A	Blue, SC1, 6mm, 23"/60cm
MMT-945A	Blue, SC1, 6mm, 32"/80cm
MMT-961A	Clear, SC1, 6mm, 18"/46cm
MMT-963A	Clear, SC1, 6mm, 23"/60cm
MMT-965A	Clear, SC1, 6mm, 32"/80cm
MMT-975A	Clear, SC1, 9mm, 32"/80cm
MMT-943T	Blue, Paradigm, 6mm, 23"/60cm
MMT-975T	Clear, Paradigm, 9mm, 32"/80cm

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5. Infusion set MiniMed Mio 30, variant:

MMT-905, MMT-905T, MMT-905A, MMT-905AT

- Insertion angle: 30 degrees
- Soft cannula outer diameter: 0.68mm
- 27 gauge introducer needle
- 13 mm cannula length
- Tubing length: 60 cm (23") and 110 cm (43").
- Tubing ID: 0.385mm, OD: 1.50mm
- Usage time: Max 72 hours

MMT-905	Mio 30, Gray, 23", 13mm cannula, Paradigm
MMT-906	Mio 30, Gray, 43", 13mm cannula, Paradigm
MMT-905T	Mio 30, Grey, 23", 13mm cannula, Paradigm
MMT-905A	Mio 30, Gray, 23", 13mm cannula, SC1
MMT-906A	Mio 30, Gray, 43", 13mm cannula, SC1
MMT-905AT	Mio 30, Grey, 23", 13mm cannula, SC1

6. Infusion set MiniMed Mio Advance, variant:

MMT-242A, MMT-242AT, MMT-213A, MMT-243A, MMT-244A

- Insertion angle: 90 degrees
- Soft cannula length: 6mm and 9mm
- Soft cannula OD ~0.68mm
- 27 gauge introducer needle
- Tubing ID: 0.385mm, OD: 1.50mm
- Tubing length: 60cm (23") and 110cm (43")
- Usage time: Max 72 hours

MMT-242A	Mio Advance, Grey, 23", 6 mm cannula, SC1
MMT-242AT	Mio Advance, Grey, 23", 6 mm cannula, SC1
MMT-213A	Mio Advance, Grey, 43", 6 mm cannula, SC1
MMT-243A	Mio Advance, Grey, 23", 9 mm cannula, SC1
MMT-244A	Mio Advance, Grey, 43", 9 mm cannula, SC1

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изделие / The list of national and international regulatory documents/standards, which matches medical device

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ISO 80369-7	Connectors with small internal diameter for liquids and gases used in medicine. Part 7. Connectors for intravascular and subcutaneous applications
ISO 80369-6	Connectors with small internal diameter for liquids and gases used in medicine. Part 6. Connectors for neuroaxial applications
ISO 2859-1	Sampling procedures for sampling by quality attributes. Part 1. Acceptable quality level (AQL) sampling plans for sequential lot control.
ISO 10555-1	Intravascular sterile single-use catheters. Part 1. General technical requirements
EN ISO 14644-1	Cleanrooms and associated controlled environments. Part 1. Classification of air cleanliness by particle concentration
EN ISO 11070	Single-use sterile intraductors, dilators and wire guides
EN 556-1	Sterilization of medical devices. Requirements for medical devices in the "sterile" category. Part 1. Requirements for sterilization of packaged medical devices.
EN ISO 11135	Sterilization of medical devices. Requirements for the development, validation and ongoing monitoring of sterilization of medical devices with ethylene oxide
EN ISO 11607-1	Packaging for medical devices subject to finishing sterilization. Part 1. Requirements for materials, sterilization barrier systems and packaging systems
EN ISO 11607-2	Packaging for medical devices subject to final sterilization. Part 2. Requirements for validation of forming, sealing and assembly processes
EN 868-5	Materials and systems for packaging of medical devices intended for sterilization. Part 5. Heat-resistant self-sealing pouches and sleeves made of porous materials and multilayer plastic film. Requirements and test methods
ISO 639-1	Codes for the representation of language names. Part 1. Two-letter code

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A ConvaTec Company

ATTESTATION STATEMENT

As a responsible member of Unomedical A/S, a Convatec company, I hereby attest that I believe the following Instructions for Use, attached herein, to be a true and accurate copies.

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





Indications for use

The MiniMed™ Mio™ Advance infusion set is indicated for subcutaneous infusion of medication administered by an external pump. The infusion set is indicated for single use.

Description

The MiniMed Mio Advance infusion set has a 90-degree soft cannula. It is delivered ready to use in a pre-loaded insertion device with automatic needle retraction. It is provided sterile and nonpyrogenic.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulin approved for subcutaneous infusion.

Contraindications

This infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous infusion. Do not use the infusion set with blood or blood products.

General instructions

- The infusion set may be worn for a maximum of three days, or per a healthcare professional's instructions.
- If the infusion set is being used for the first time, do the first set-up in the presence of a healthcare professional.

Warnings

- Do not use the infusion set if the package is opened or damaged. An opened or damaged package may be contaminated. This can cause

an infection.

- Do not use the infusion set if the disconnect cover is removed. Use a new infusion set instead. One purpose of the disconnect cover is as a safety feature against accidental firing of the insertion device which could cause injury.
- Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours afterwards.
- Do not put alcohol, disinfectants, perfumes, deodorants, cosmetics, or other substances with solvents on the infusion set. These substances may damage the integrity of the infusion set.
- Never point the insertion device toward any body part where insertion is not desired.
- Always rotate the insertion site when changing the infusion set. Reusing an insertion site too often may cause scarring and unpredictable insulin delivery. Review the pump user guide on how to rotate insertion sites.
- Check the insertion site often through the clear window. Improper insertion and maintenance of the insertion site can cause inaccurate insulin delivery, infection, or site irritation. Replace the infusion set at a new site if the soft cannula is not properly inserted.
- Check the cannula housing and tubing for blood. Blood may cause insufficient insulin delivery. This may result in high blood glucose. Replace the infusion set at a new site if there is blood.
- Fill the tubing completely with insulin before insertion. Do not leave air in the tubing. An uncontrolled amount of insulin may be delivered if air is present.
- Check the tubing for clogs or leaks if a blood glucose reading is high. Clogs or leaks may

restrict insulin delivery and result in high blood glucose. Replace the infusion set even if a clog or leak is suspected, but not found.

- Do not attempt to clear air or an occlusion in the tubing while it is connected to the body. An uncontrolled amount of insulin may be delivered. This can cause high or low blood glucose. Disconnect the tubing before adjusting it.
- Do not reuse the infusion set. Reuse of the infusion set may cause damage to the cannula or needle and may lead to infection, site irritation, or unpredictable insulin delivery.
- If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia. If this occurs, start over with a new reservoir and infusion set.
- For no delivery alarms, refer to the pump user guide.

Precautions

- Do not use another type of infusion set without consulting a healthcare professional for correct handling. A healthcare professional should always be consulted when choosing an infusion set.
- Choose insertion sites as recommended by a healthcare professional. The choice of site depends on treatment and patient-specific factors such as body composition and physical activity level.
- Check that the insertion site is free of skin irritation, such as redness, scar tissue, or bleeding. Do not insert the infusion set into

muscle or over bone. Doing so can cause pain or may damage the infusion set. If the infusion set is improperly inserted, replace the infusion set and insert at a new site.

- If needed, remove body hair around the insertion site to ensure that the adhesive sticks to the skin.
- Do not use an insertion site located under a belt or waistband, or where the site is constrained by clothing or accessories. Insertion at these locations can cause disconnection of the infusion set and interruption of insulin delivery, and lead to hyperglycemia.
- Do not reposition the infusion set on the body after the adhesive is placed on the skin. Repositioning the infusion set may damage the adhesive. Replace the infusion set if the adhesive is damaged.
- Always check blood glucose one to three hours after inserting a new infusion set. This is to confirm accurate insulin delivery. High or low blood glucose can result if insulin delivery is not accurate.
- Check blood glucose several times throughout the day or as recommended by a healthcare professional.
- Consult a healthcare professional for how to correct for missed insulin and for how much time the pump can remain disconnected.

Storage and disposal

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the

manufacturer's instructions.

- Dispose of the insertion device in a proper sharps container according to local laws.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

Serious incident

- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and to your national authority.

Warranty


For product warranty information, please contact a local Medtronic support representative, or visit: www.medtronicdiabetes.com/warranty

Additional information for the healthcare professional

- Consider treatment and patient-specific factors when choosing an insertion site.
- The length of the soft cannula should be based on treatment and patient factors such as physiology and activity level. If the soft cannula is too long, insertion pain or insertion into the bone or muscle may occur. If the soft cannula is too short, leakage or irritation at the insertion site may occur. In both cases, the infusion must be stopped and a new soft cannula with a different length must be inserted at a new insertion site.
- Inform the patient that the selection of an insertion site is dependent on cannula length.

©2019 Medtronic MiniMed, Inc. All rights reserved. MiniMed and Mio are trademarks of Medtronic MiniMed, Inc.

Components

- | | |
|--|--|
| A. Top button | H. Paper backing |
| B. Disconnect cover | I. Adhesive |
| C. Insertion device | J. MiniMed or Luer-lock tubing connector |
| D. Raised arrows indicating tubing direction | K. Tubing |
| E. Clear Window | L. Site Connector |
| F. Cannula housing | M. White cap for site connector |
| G. Soft cannula | |

Instructions for use

- Read all instructions carefully before using the MiniMed Mio Advance infusion set.
- Consult the pump user guide for important information about pump therapy prior to connecting the infusion set. This information

includes connection and filling procedures, possible errors, and potential risks related to pump therapy.

- Follow good hygienic procedures. If this infusion set is being used for the first time, perform the first set-up in the presence of a healthcare professional.

Inserting

Unfold the front cover of this booklet to view images that are keyed to the steps below.

1. Wash hands with soap and water.
2. Select a recommended insertion site (shown in gray) as indicated by a healthcare professional.
3. Clean the insertion site with a disinfectant as directed by a healthcare professional. Allow to air dry before inserting the infusion set. Remove hair around the insertion site to ensure the adhesive sticks to the skin if needed.
4. Open the package and remove the paper from the tubing.
5. Place the tubing connector on the top of a newly filled reservoir. If a MiniMed tubing connector is in use, twist the tubing connector clockwise until it locks in place. If a Luer-lock tubing connector is in use, twist the tubing connector until it is tightly secured.

For the MiniMed tubing connector, make sure there is no liquid inside the tubing connector or on top of the reservoir. Liquid can block vents and cause inaccurate insulin flow. This can result in high or low blood glucose. If there is

liquid in the vents, use a new infusion set and a new reservoir.

6. Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Keep the white cap for later use.
7. Put the reservoir into the pump, then completely fill the tubing with insulin as instructed by the pump user guide. Do not leave any air in the tubing. The tubing is completely filled when drops of insulin are seen coming out of the site connector.
8. Remove the paper backing from the adhesive. Be careful not to touch the adhesive.
9. Remove the disconnect cover from the insertion device by (a) gently squeezing the sides of the disconnect cover then (b) pulling it away from the insertion device. Keep the disconnect cover for later use. The raised arrows indicate the tubing direction.
- 10a. Stretch the skin until smooth. Then press the insertion device against the skin.
- 10b. Press the top button down completely to insert the infusion set.
11. Gently and carefully remove the insertion device.
12. Gently massage the adhesive onto the skin with one finger. Replace the infusion set if the adhesive does not stick to the skin.
13. Gently hold the cannula housing steady with one finger. Then push the site connector straight into the cannula housing until it clicks.

14. Fill the soft cannula with insulin:

6 mm: 0.6 units (0.006 ml)

9 mm: 0.6 units (0.006 ml)

Dispose of the insertion device in an appropriate sharps container according to local laws.

Keep the disconnect cover and the white cap for use when the infusion set is disconnected.

Disconnecting

1. Wash hands with soap and water.
2. Gently hold the cannula housing steady with one finger. Then (a) squeeze the sides of the site connector and (b) pull the site connector out of the cannula housing.
3. Put the white cap on the site connector.
4. Put the disconnect cover on the cannula housing. Then push the cover into the cannula housing until it clicks.

Reconnecting

1. Wash hands with soap and water.
2. Use one finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the disconnect cover and (b) pull the disconnect cover out of the cannula housing.
3. Remove the white cap from the site connector by (a) squeezing the sides of the site connector then (b) pulling the white cap off. Ensure that there is no air in the tubing. If there are air bubbles in the tubing: Fill the tubing with insulin as instructed in the pump user guide. Do not leave any air in the tubing.
4. Use one finger to gently hold the cannula

housing steady. Then, push the site connector into the cannula housing until it clicks.

Removing

1. Wash hands with soap and water.
2. Use one finger to gently hold the cannula housing steady. Then (a) squeeze the sides of the site connector and (b) pull the site connector out of the cannula housing.
3. Disconnect the tubing from the pump and remove the tubing connector from the reservoir by turning the connector counter-clockwise. For specific instructions to disconnect tubing from the pump, refer to the pump user guide.
4. Carefully lift the adhesive around the cannula housing. Then, pull the soft cannula out of the skin.

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Contacts:**Africa:**

Medtronic South Africa and Southern
Africa

Office Reception Tel:
+27(0) 11 260 9300

Diabetes: 24/7 Helpline:
0800 633 7867

Sub-Sahara 24/7 Helpline:
+27(0) 11 260 9490

Albania:

Net Electronics Albania
Tel: +355 697070121

Argentina:

Corpomedica S.A.
Tel: +(11) 4 814 1333
Medtronic Directo 24/7:
+0800 333 0752

Armenia:

Exiol LLC
Tel: +374 98 92 00 11
or +374 94 38 38 52

Australia:

Medtronic Australasia Pty. Ltd.
Tel: 1800 668 670

Bangladesh:

Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-8903995417
or (+880)-1714217131

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





















Indications for use

The MiniMed Mio infusion set is indicated for the subcutaneous infusion of insulin from an infusion pump. The infusion set is indicated for single use.

Intended use

The infusion set is indicated for subcutaneous infusion of insulin in the treatment of diabetes mellitus. The infusion set has been tested to be compatible with insulins approved for subcutaneous infusion.

Description

The MiniMed Mio is an infusion device which combines an infusion set (catheter) with an aid for insertion (serter). As the infusion set is an integrated part of theserter, the infusion set is delivered assembled and ready for use as a single, sterile unit.

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Warnings

- Read all instructions carefully before using the infusion set. Failure to follow instructions may result in pain or injury.
- The infusion set is sterile and non-pyrogenic unless the package has been opened or damaged. Do not use the infusion set if the package has been opened or damaged.

- For Medtronic pump users, do not use the infusion set if the tubing connector needle has been damaged.
- Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged.
- If using this infusion set for the first time, do the first set-up in the presence of a healthcare professional.
- Inaccurate medication delivery, site irritation or infection may result from improper insertion and maintenance of the infusion site.
- Do not re-insert the introducer needle into the infusion set. This can cause tearing or damage to the soft cannula and can lead to unpredictable medication flow.
- Reuse of the infusion set can cause site irritation or infection and can lead to inaccurate medication delivery.
- Always remove the needle guard before inserting the infusion set.
- Before insertion, clean the insertion site with a disinfectant as indicated by a healthcare professional.
- For the MiniMed tubing connector: If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

Precautions

- Never try to fill or free clogged tubing while the infusion set is inserted. This can result in unpredictable medication flow.
- When unwinding the tubing, release it with caution, because a hard pull on the tubing can result in damage to the infusion set or introducer needle. Ensure that the infusion set is properly in place when the tubing is fully released.
- Do not leave air in the infusion set. Make sure to fill the tubing completely. Consult the pump manufacturer's instructions on how to fill the tubing.
- The infusion set is a single use item, which must be disposed of after use. Do not clean or re-sterilize.
- Do not put disinfectants, perfumes, or deodorants on the infusion set as these may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose, or if the soft cannula is displaced from its original site.
- Check the infusion site frequently to ensure that the soft cannula remains firmly in place. If the infusion set is not in place, replace it with a new one at a new site. Since the cannula is soft, it will not cause any pain if it slips out, which may happen without notice. The soft cannula must always be completely inserted to receive the full amount of medication.
- Replace the infusion set every two to three days, or per a healthcare professional's instructions.

- If the infusion site becomes irritated or infected, replace the set, and use a new site until the first site has healed.
- Use aseptic techniques when temporarily disconnecting the infusion set. Consult a healthcare professional for how to compensate for missed medication while disconnected.
- Carefully monitor blood glucose levels when disconnected, and after reconnecting the infusion set.
- Never point a loaded insertion device towards any body part, where insertion is not desired.
- Check blood glucose levels one to three hours after inserting the infusion set. Measure blood glucose on a regular basis. Discuss this with a healthcare professional.
- Do not change the infusion set just before bedtime unless blood glucose can be checked one to three hours after insertion.
- If blood glucose becomes inexplicably high, or an occlusion alarm occurs, check for clogs or leaks. If in doubt, change the infusion set, because the soft cannula may be dislodged or partially clogged. Discuss a plan to lower blood glucose levels with a healthcare professional. Test blood glucose to ensure the problem has been corrected.
- If, during the use of this device or as a result of its use a serious incident has occurred, please report it to the manufacturer and your national authority.

Storage and disposal

- Store infusion sets in a cool, dry location at room temperature. Do not store infusion sets in direct sunlight, or in high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid risk of needle stick injury.
- Dispose of a used infusion set according to local regulations for biohazardous waste.


Warranty

- For product warranty information, please contact a local Medtronic support representative or visit: www.medtronicdiabetes.com/warranty.

ENGLISH

Preparation and insertion

This product is an infusion device which combines an infusion set (catheter) with an aid for insertion (serter). This single-use, sterile product is delivered assembled and ready for use. Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged. Carefully examine the instructions for use before inserting this product. Refer to pump and reservoir user guides for installation and filling. Please see the corresponding illustrations on pages 2 and 3.

Contents

- A. Lid
- B. Needle guard
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Site connector
- G. Tubing

- H. Introducer needle
- I. Insertion device
- J. Sterile paper

Important note!

- K. Three raised dots – press here when lifting the lid off (see Figure 6)
- L. Lined indentations – press here when loading (see Figure 10)
- M. Round Indentations – press both sides when releasing (see Figure 12)
- N. Disconnect cover (see Figure B)

Instructions

1. Wash your hands.
2. Suggested infusion sites. Choose infusion sites recommended by a healthcare professional. Do not use an area right next to your previous insertion site.
3. Clean the insertion area with a disinfectant as directed by your healthcare professional. Be sure the area is air dry before inserting the infusion set.
4. Pull the red tape to remove the protective seal.
5. Remove the sterile paper.
6. Hold the bottom of the insertion device with one hand, while firmly pressing the three raised dots on each side of the lid with the other hand, and lift up the lid.
CAUTION: Be careful not to bend or touch the introducer needle when preparing theserter and infusion set.

7. **Unwinding the tubing:** Carefully twist the beginning of the tubing out of the slot. Now unwind the tubing from the insertion device by pulling it carefully upwards.

Caution: Do not touch the introducer needle when unwinding the tubing. Be careful not to pull too hard when the tubing is fully released. You could accidentally remove the infusion set from the introducer needle.

IMPORTANT NOTE! If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead prior to insertion, pain or minor injury may occur. To avoid this, please ensure that the following steps are taken

gently so the infusion set is not loosened / removed by accident.

8. Connect the tubing to a filled reservoir. Fill the infusion set. Ensure that there are no air bubbles in the tubing. Consult the pump manufacturer's instructions for use.

Caution: When filling the infusion set, hold the infusion set with the needle pointing down to ensure that medication does not come in contact with the adhesive backing paper.

9. Pull up gently to remove the backing paper.
10. **To prepare theserter:** Place fingers on lined indentations. Pull up the spring until you hear a "click".
11. Carefully remove the needle guard by gently twisting and then pulling it off. Make sure that the soft cannula does not exceed the introducer needle.
Caution: Before positioning the infusion set over an insertion site, place the tubing in the slot on the side of theserter to ensure that the tubing is not caught underneath the infusion set during insertion.
12. **To insert the infusion set:** Position the infusion set over the insertion site. Simultaneously press the round indentations on both sides of the insertion device to insert the infusion set.
13. Push gently on the spring to secure the adhesive onto the skin.
14. Remove theserter and introducer needle by grasping the center of theserter and pulling it gently back. Gently massage the adhesive tape

onto the skin. Fill the empty space in the cannula with the following amount of U-100 insulin:

U-100 insulin:

6 mm: 0.3 units (0.003 ml)

9 mm: 0.5 units (0.005 ml)

15. Push the lid back in place until you hear a "click."

To disconnect the infusion set

The infusion set allows you to temporarily disconnect from your pump without having to change the infusion set.

- A. Hold the adhesive tape in place by placing a finger just in front of the cannula housing, and gently squeeze the sides of the site connector. Pull the site connector straight out from the cannula housing.
- B. Insert the disconnect cover into the cannula housing until you hear a "click".

To reconnect the infusion set

- C. If necessary, fill the infusion set until medication dispenses from the needle.
Caution: When filling the infusion set, hold the needle pointing down.
- D. Remove the disconnect cover from the cannula housing by pressing the grips. Place a finger on the adhesive tape just in front of the cannula housing.
- E. Push the connector needle straight in until you hear a "click"

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Medtronic Australasia Pty. Ltd.
Tel: 1800 668 670

Bangladesh:

Sonargaon Healthcare Pvt Ltd.
Mobile: (+91)-9903995417
or (+880)-1714217131

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USA:

Medtronic Diabetes Global
Headquarters
24-Hour Technical Support:
+1-800-646-4633
To order supplies: +1-800-843-6687

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

Indications for use

The MiniMed Mio 30 infusion set is indicated for the subcutaneous infusion of insulin from an infusion pump. The infusion set is indicated for single use.

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Description

MiniMed Mio 30 is an angled infusion set with an integrated insertion device. The built-in legs support a 30° insertion angle. MiniMed Mio 30 is delivered ready for use. Both MiniMed Mio 30 and the separate tubing unit are sterile.

Contraindications

The infusion set is indicated for subcutaneous use only. Do not use the infusion set for intravenous (IV) infusion. Do not use the infusion set for infusion of blood or blood products.

Warnings

- Read all instructions for use before inserting the infusion set. Failure to follow instructions may result in pain or injury
- The infusion set is only sterile and non-pyrogenic if the package is undamaged and unopened. Do

not use the infusion set if the package is open or damaged. Ensure sterility by checking that the sterile paper and tamper-proof seal are not damaged.

- If using the infusion set for the first time, do so in the presence of a healthcare professional.
- Inaccurate medication delivery, site irritation, or infection may result from improper insertion and maintenance of the infusion site.
- Do not reinsert the introducer needle into the infusion set. Reinsertion may cause damage or tearing of the soft cannula, which may lead to unpredictable medication flow.
- Do not reuse the infusion set. Reuse of the infusion set may cause site irritation or infection and may lead to inaccurate medication delivery.
- Carefully remove the needle guard before inserting the infusion set. Do not bend the needle prior to insertion. Do not use the infusion set if the needle is bent or has been damaged.
- If insulin or any liquid gets inside the tubing connector, it can temporarily block the vents that allow the pump to properly fill the infusion set. **This may result in the delivery of too little or too much insulin, which can cause hyperglycemia or hypoglycemia.** If this occurs, start over with a new reservoir and infusion set.

Precautions

- Never try to fill or free a clogged tubing while the infusion set is inserted. This may result in unpredictable medication flow.
- Do not leave air in the infusion set or the tubing. Make sure to fill the tubing completely. For filling instructions, consult the pump manufacturer's instructions for use.
- The infusion set is a single-use device which must be disposed of after use. Do not clean or resterilize.
- Check the infusion set frequently to ensure that the soft cannula remains firmly in place. Since the cannula is soft, it will not cause any pain if it slips out, and this can happen without notice. Replace with a new infusion set at a new site if the infusion set is not in place. The soft cannula must always be completely inserted to receive the full amount of medication.
- Do not put disinfectants, perfumes, alcohol, or deodorants on the infusion set, as these may impact the integrity of the infusion set.
- Replace the infusion set if the adhesive tape becomes loose, or if the soft cannula is displaced from its original site. If the tape is curled, do not attempt to straighten it. Stretching or pulling on the tape may affect the cannula under the skin.

- Replace the infusion set every two to three days, or as directed by a healthcare professional.
- If the infusion site becomes irritated or infected, replace the infusion set, and use a new site until the previous site has healed.
- Use aseptic techniques when temporarily disconnecting the infusion set. Consult a healthcare professional for how to compensate for missed medication while the infusion set is disconnected.
- To avoid accidental needle stick injury, never point a loaded insertion device toward any body part where insertion is not desired.
- Check blood glucose one to three hours after insertion. Measure blood glucose regularly. Consult a healthcare professional for more information.
- Do not change the infusion set just before bedtime, unless blood glucose can be checked one to three hours after insertion.
- If blood glucose becomes inexplicably high, or an occlusion alarm occurs, check the infusion set for clogs or leaks. If in doubt, change the infusion set, because the soft cannula may be dislodged. Discuss a plan with a healthcare professional to lower blood glucose. Test blood glucose to make sure the problem is corrected.

- Carefully monitor blood glucose levels when disconnected and after reconnecting the infusion set.
- Disconnect the tubing with caution by pressing the grips. A hard pull on the tubing can result in damage to the infusion set. Ensure that the infusion set is properly in place when the tubing is fully disconnected.

Serious incident

- If, during the use of this device, or as a result of its use a serious incident occurs, report the incident to the manufacturer and to the national authority.

Storage and disposal

- Store the infusion set in a cool, dry location at room temperature. Do not store the infusion set in direct sunlight or high humidity.
- Store and handle insulin according to the manufacturer's instructions.
- Put the lid back on the insertion device and dispose in a proper sharps container according to local laws, to avoid the risk of needle stick injury.
- Dispose of a used infusion set according to local regulations for biohazardous waste.

Warranty

For product warranty information, contact a local Medtronic support representative, or visit: www.medtronicdiabetes.com/warranty.

© 2020 Medtronic. Medtronic, Medtronic with rising man logo, and Medtronic logo are trademarks of Medtronic. Third-party trademarks (TM) belong to their respective owners. The following list includes trademarks or registered trademarks of a Medtronic entity in the United States and/or in other countries.


MiniMed™ Mio™ 30

MiniMed Mio 30 instructions for use

Read all instructions carefully before using MiniMed Mio 30. Consult the pump user guide before connecting MiniMed Mio 30 to the pump. Follow proper hygiene procedures. If using MiniMed Mio 30 for the first time, do the first set-up in the presence of a healthcare professional.

Contents

- A. Lid
- B. Needle guard
- C. Soft cannula
- D. Adhesive tape
- E. Cannula housing
- F. Site connector
- G. Tubing connector
- H. Introducer needle
- I. Insertion device
- J. Sterile paper
- K. Finger grips – for holding MiniMed Mio 30

- L. Release buttons – press here when releasing MiniMed Mio 30
- M. Disconnect cover
- N. Protective cap for site connection

- O. Legs
- P. Spring with finger loop (inside)

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CAUTION: Check the tamper-proof seal and sterile paper. Do not use if the seal or paper has been broken.

1. Wash your hands before inserting the infusion set.
2. Select a recommended insertion site, shown in grey, and as indicated by a healthcare professional. Consult a healthcare professional for insertion site rotation.
3. Clean the area of insertion with a disinfectant as directed by a healthcare professional. Allow skin to air dry for at least 30 seconds before inserting the infusion set. If needed, remove hair around the insertion site to prevent improper attachment of the tape to the skin.
4. Remove the paper from the tubing. Connect the tubing to a filled reservoir.

IMPORTANT: Make sure there is no liquid inside the tubing connector or on top of the reservoir. If this occurs, start over with a new reservoir and infusion set.

5. Remove the protective cap of the site connector. Fill the tubing until insulin drips from the tip. Ensure there are no air bubbles in the tubing. Consult the pump manufacturer's instructions for use for proper filling procedures.

6. Pull the red tape to remove the seal.
7. Pull off the sterile paper. Do not touch the spring at this stage.
8. Hold the MiniMed Mio 30 with one hand and pull the lid off with the other.
CAUTION: Be careful not to bend or touch the MiniMed Mio 30 introducer needle. If the infusion set is not securely placed in the insertion device with the needle pointing straight ahead, pain or minor injury may occur during insertion.
9. Place your fingers on the finger grips on each side. With your other hand, use your finger to gently pull back the spring until you hear a "click". The needle guard will lift up or peel off. Do not touch or bend the introducer needle. Check that the soft cannula does not exceed the introducer needle.
10. Be sure the adhesive tape is not stuck on the introducer needle.
11. **INSERTION:** Place your index finger on the top release button and your thumb on the bottom release button. Position MiniMed Mio 30 with its legs resting flat on the skin to ensure an insertion angle of 30°. Make sure to maintain a 30° angle while inserting and gently press the release buttons once to insert MiniMed Mio 30.
12. Make sure that the cannula stays in place. Gently use two fingers to hold down the tape while

removing the insertion device and introducer needle. Gently pull the insertion device straight back.

13. Secure the cannula by gently pressing down with one finger on the clear window of the tape. Carefully remove the protective liner under the cannula housing.
14. Remove the backing paper under the front end. Gently massage the adhesive tape onto the skin.
15. Put the lid back onto the MiniMed Mio 30 for disposal.
16. Gently place a finger on the cannula housing to hold it and insert the site connector into the infusion set using your other hand. Make sure the connector locks with a "click". Fill the empty space in the cannula with 0.7 units (0.007 ml) of U-100 insulin.

Disconnecting

MiniMed Mio 30 allows you to temporarily disconnect your pump without having to change the infusion set.

- A. Gently place a finger in front of the cannula housing and squeeze both sides of the site connector. Pull the connector straight out.
- B. Insert the protective cap onto the site connector until you hear a "click".
- C. Insert the disconnect cover into the cannula housing until you hear a "click".

Reconnecting

- D. Remove the disconnect cover from the cannula housing.
- E. Remove the protective cap from the site connector. If necessary, fill the infusion set until insulin dispenses from the needle.
CAUTION: When filling, hold MiniMed Mio 30 with the needle pointing down.
- F. Gently place a finger in front of the cannula housing. Push the site connector into the cannula housing until you hear a "click".

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Contacts:

Africa:

Medtronic South Africa and Southern Africa

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


























This is to certify that Mr Michael Christian Remy Jensen today in my presence at the Notarial Office approved and signed the above document. No conspicuous corrections or addenda were found in the document.

Mr Michael Christian Remy Jensen has proved his identity by presenting Danish Passport No 209696572.

District Court of Lyngby, Denmark, 21 April 2023

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