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Description: Avitene™ UltraFoam™ is an absorbable hemostatic sponge prepared as a sterile, porous, pliable, water insoluble partial hydrochloric acid salt of purified bovine corium collagen. In its manufacture, swelling of the native collagen fibrils is controlled by ethyl alcohol to permit noncovalent attachment of hydrochloric acid to amine groups on the collagen molecule and preservation of the essential morphology of native collagen molecules. The characteristics of collagen which are essential to its effect on the blood coagulation mechanisms are preserved, although dry heat sterilization causes some cross-linking which is evidenced by reduction of hydrating properties, and a decrease of molecular weight which implies some degradation of collagen molecules.

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Actions: Avitene™ UltraFoam™ sponge, in contact with a bleeding surface, attracts platelets which then undergo the release phenomenon to trigger aggregation of platelets into thrombi in the interstices of the fibrous mass. It is able to absorb and hold within its interstices many times its weight of blood and other fluids. When implanted in animal tissues, it is absorbed in 90 days.

Indications: Avitene™ (MCH) is used in surgical procedures as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

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Warnings:

- Avitene™ UltraFoam™ sponge is inactivated by autoclaving.
- Ethylene oxide reacts with bound hydrochloric acid to form ethylene chlorohydrin.
- This device has been designed for single use only. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient.

Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may

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lead to injury, illness or death of the patient or end user. Opened, unused product should be discarded.

- As with any foreign substance, use of Avitene™ UltraFoam™ sponge in contaminated wounds may enhance infection.
 - Avitene™ UltraFoam™ sponge should not be used in instances of pumping arterial hemorrhage.
 - It should not be used where blood or other fluids have pooled, or in cases where the point of hemorrhage is submerged as it may mask an underlying source of bleeding, resulting in hematoma.
 - Avitene™ UltraFoam™ sponge will not act as a tampon or plug in a bleeding site, nor will it close off an area of blood collecting behind a tampon.
- Avitene™ UltraFoam™ sponge is not intended to treat systemic coagulation disorders.
- Not for injection, intraocular or intravascular use.

Precautions: • Avitene™ UltraFoam™ sponge, like other collagen hemostats, should not be left in infected sites or on infected surfaces.

- As with other hemostatic agents, it is not recommended for use in persons with known sensitivity to material of bovine origin.
- Any excess Avitene™ UltraFoam™ not removed at the time of surgery may either present itself as a (recurring) mass or a (space occupying) lesion or it may lead to a foreign body reaction that may present with or without clinical signs and symptoms as a recurring mass or lesion or postoperative abscess formation upon imaging. Imaging may initially not be capable of distinguishing the difference. Removal of excess material, ideally performed upon conclusion of the initial procedure, typically resolves all signs and symptoms.
- Microfibrillar Collagen Hemostat (MCH) contains a low, but detectable, level of intercalated bovine serum protein which reacts immunologically as does beef serum albumin. Increases in anti-BSA titer have been observed following treatment with MCH. About two-thirds of individuals exhibit antibody titers because of ingestion of food products of bovine origin. Intradermal skin tests have occasionally shown a weak positive reaction to BSA or MCH but these have not been correlated with IgG titers to BSA. Tests have failed to demonstrate clinically significant elicitation of antibodies of the IgG class against BSA following MCH therapy.
- When placed into cavities or closed spaces, care should be exercised to avoid overpacking Avitene™ UltraFoam™ sponge as it may press against neighboring structures.
- The safety of this product has not been established in children or in pregnant women; therefore, Avitene™ UltraFoam™ sponge should only be used after an evaluation of the relative benefits and risks clearly warrant its use.

Adverse Reactions:

- The most serious adverse reactions reported, which may be related to the use of other collagen products, are potentiation of infection including abscess formation, hematoma, wound dehiscence and mediastinitis.
- Other reported adverse reactions possibly related are adhesion formation, allergic reaction, foreign body reaction and subgaleal seroma (report of a single case), and increased incidence of alveolalgia when used for packing of dental extraction sockets.

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Product Form	Product Code	Characteristics
Avitene™ Microfibrillar Collagen Hemostat UltraFoam™	1050020	Size: 2cm x 6.25cm x 7mm Mass: not more than 0,2 g Percent solids: 1.10 -1.64; nominal = 1.37%

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		Cast Ultrafoam: 0.85% solids (9.5 grams of Avitene per liter of USP water) Modulus (kPa): not more than 592,949 Acute Max. Load: not less than 0,36 N (36,29 g)
	1050030	Size: 8cm x 6.25cm x 1cm Mass: not more than 0,5 g Percent solids: 1.10 -1.64; nominal = 1.37% Cast Ultrafoam: 0.85% solids (9.5 grams of Avitene per liter of USP water) Modulus (kPa): not more than 592,949 Acute Max. Load: not less than 0,36 N (36,29 g).
	1050040	Size: 8cm x 12.5cm x 1cm Mass: not more than 0,5 g Percent solids: 1.10 -1.64; nominal = 1.37% Cast Ultrafoam: 0.85% solids (9.5 grams of Avitene per liter of USP water) Modulus (kPa): not more than 592,949 Acute Max. Load: not less than 0,36 N (36,29 g)
	1050050	Size: 8cm x 12.5cm x 3mm Mass: not more than 0,5 g Percent solids: 1.10 -1.64; nominal = 1.37% Cast Ultrafoam: 0.85% solids (9.5 grams of Avitene per liter of USP water) Modulus (kPa): not more than 592,949 Acute Max. Load: not less than 0,36 N (36,29 g).

Tolerance is $\pm 5\%$ if unless otherwise specified

Package contents:

Avitene™ UltraFoam™ sponge is supplied sterile and non-pyrogenic. The sterility of this product is guaranteed unless the foil pouch is opened or damaged.

Product Form	Product Codes	Units/case
Avitene™ Microfibrillar Collagen Hemostat UltraFoam™	1050020	12 unit per package
	1050030	6 unit per package
	1050040	6 unit per package
	1050050	6 unit per package

Dosage And Administration: Only the amount of material necessary to control bleeding should be used. Avitene™ UltraFoam™ sponge should be cut to the desired size and applied directly to the source of bleeding. Hold Avitene™ UltraFoam™ sponge in place with moderate pressure until hemostasis results. The period of time to hold pressure will vary with the force and severity of bleeding. Avitene™ UltraFoam™ sponge may be left in place at the bleeding site when necessary. Prior to removal, moisten Avitene™ UltraFoam™ sponge with saline to avoid dislodging the clot.

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Potential consumer: Medical device should only be used by qualified medical personnel with experience in the use of similar products or by medical personnel who have completed training in working with similar products under the guidance of an experienced specialist with experience in the use of similar products. Avitene MCH is a prescription device and as such should be used

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by a qualified medical professional.

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Information about materials of animal/human origin: Avitene™ Microfibrillar Collagen Hemostat does incorporate or utilize any substance or material that is derived from human or animal origin as defined in Article 1.5(f)(g) of Directive 93/42/EEC.

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Shelf life/Service life: The Avitene™ UltraFoam™ is labeled with a 38 month shelf life. Duration of absorption in human body - 8-10 weeks.

Disposal: After use, the applicators are potentially biohazardous products due to exposure to blood and body fluids (class B medical devices under medical wastes classification). Dispose of any potentially biohazardous products in accordance with accepted medical practice and

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
applicable local, state and federal regulations and the institutions written policy for the disposal of potentially biohazardous material.

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Information about manufacturer:**Legal Manufacturer**

Davol, Inc. Subsidiary of C. R. Bard, Inc.
100 Crossings Blvd.
Warwick RI, 02886
USA

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State of Rhode Island
County of Providence

On this 23rd day of March, 2020, before me, the undersigned notary public, personally appeared Penny Cardoso, and proved to me through satisfactory evidence of identification to be the person whose name is signed on the preceding or attached document, and acknowledged that they signed it voluntarily for its stated purpose.

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Description: AVITENE™ (Microfibrillar Collagen Hemostat, or MCH), is an absorbable topical hemostatic agent prepared as a dry, sterile, fibrous, water insoluble partial hydrochloric acid salt of purified bovine corium collagen. It is prepared in a loose fibrous form and in a compacted “non-woven” web form. In its manufacture, swelling of the native collagen fibrils is controlled by ethyl alcohol to permit non-covalent attachment of hydrochloric acid to amine groups on the collagen molecule and preservation of the essential morphology of native collagen molecules. Dry heat sterilization causes some cross-linking which is evidenced by reduction of hydrating properties, and a decrease of molecular weight which implies some degradation of collagen molecules. However, the characteristics of collagen, which are essential to its effect on the blood coagulation mechanisms, are preserved.

Purpose of the medical device: Avitene™ Microfibrillar Collagen Hemostat (MCH) is indicated in surgical procedures as an adjunct to hemostasis when control of bleeding by ligature or conventional procedures is ineffective or impractical.

Actions: Avitene™ (MCH), in contact with a bleeding surface, attracts platelets which adhere to the fibrils and undergo the release phenomenon to trigger aggregation of platelets into thrombi in the interstices of the fibrous mass. The effect on platelet adhesion and aggregation is not inhibited by heparin in vitro. It has been found effective in heparinized dogs and in eight of nine fully heparinized human subjects. Platelets of patients with clinical thrombasthenia do not adhere to MCH in vitro. However, in clinical trials it was effective in 50 of 68 patients receiving aspirin. MCH cannot control bleeding due to systemic coagulation disorders. Appropriate therapy to correct the underlying coagulopathy should be instituted prior to use of the product. MCH is tenaciously adherent to surfaces wet with blood but excess material not involved in the hemostatic clot may be removed by teasing or irrigation, usually without causing rebleeding. In animal and human studies, it has been shown to stimulate a mild, chronic cellular inflammatory response. When implanted in animal tissues, it is absorbed in less than 84 days. In human studies of hemostasis in osteotomy cuts, it has been shown not to interfere with bone regeneration or healing. In animal studies, it has been demonstrated that MCH does not predispose to stenosis at vascular anastomotic sites. These findings have not been confirmed in human use. Studies have been performed using MCH (fibrous form) in experimental wounds contaminated with hemolytic *Staphylococcus aureus*. The presence of MCH does not enhance or initiate *Staphylococcus* wound infections to a greater or lesser extent than control agents employed for the same purpose.

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

Indications: Avitene™ (MCH) is used in surgical procedures as an adjunct to hemostasis when control of bleeding by ligation or conventional procedures is ineffective or impractical.

Contraindications: Avitene™ (MCH) should not be used in the closure of skin incisions as it may interfere with the healing of the skin edges. This is due to simple mechanical interposition of dry collagen and not to any intrinsic interference with wound healing. By filling porosities of cancellous bone, MCH may significantly reduce the bond strength of methylmethacrylate adhesives. MCH should not, therefore, be employed on bone surfaces to which prosthetic materials are to be attached with methylmethacrylate adhesives.

Warnings: Avitene™ (MCH) is inactivated by autoclaving. Ethylene oxide reacts with bound hydrochloric acid to form ethylene chlorohydrin. This device has been designed for single use only. Reuse, reprocessing, resterilization or repackaging may compromise the structural integrity and/or essential material and design characteristics that are critical to the overall performance of the device and may lead to device failure which may result in injury to the patient. Reuse, reprocessing, resterilization or repackaging may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient or end user. MCH is not for injection or intraocular use. Moistening MCH or wetting with saline or thrombin impairs its hemostatic efficacy. It should be used dry. Discard any unused portion. As with any foreign substance, use in contaminated wounds may enhance infection.

Warning: Not for injection or intravascular use.

Precautions: Only that amount of Avitene™ (MCH) necessary to produce hemostasis should be used. After several minutes, excess material should be removed; this is usually possible without the reinitiation of active bleeding. Any excess Avitene™ (MCH) not removed at the time of surgery may either present itself as a (recurring) mass or a (space occupying) lesion or it may lead to a foreign body reaction that may present with or without clinical signs and symptoms as a recurring mass or lesion or postoperative abscess formation upon imaging. Imaging may initially not be capable of distinguishing the difference. Removal of excess material, ideally performed upon conclusion of the initial procedure, typically resolves all signs and symptoms. Failure to remove excess MCH may result in bowel adhesion or mechanical pressure sufficient to compromise the ureter. In otolaryngological surgery, precautions against aspiration should include removal of all excess dry material and thorough irrigation of the pharynx. MCH contains a low, but detectable, level of intercalated bovine serum protein which reacts immunologically as does beef serum albumin. Increases in anti-BSA titer have been observed en following treatment with MCH. About two-thirds of individuals exhibit antibody titers because of ingestion of food products of bovine origin. Intradermal skin tests have occasionally shown a weak positive reaction to BSA or MCH but these have not been correlated with IgG titers to BSA. Tests have failed to demonstrate clinically significant elicitation of antibodies of the IgG class against BSA following MCH therapy. Care should be exercised to avoid spillage on nonbleeding surfaces particularly in abdominal or thoracic viscera. Avitene™ (MCH) should not be used in conjunction with autologous blood salvage circuits, as Avitene™ may pass through the filters of such systems. It has been suggested that fragments of MCH may pass through filters of blood scavenging systems, therefore the reintroduction of blood from operative sites treated with MCH should be avoided. Teratology studies in rats and rabbits have revealed no harm to the animal fetus. There are no well-controlled studies in pregnant women, therefore, MCH should be used in pregnant women only when clearly needed. Avitene™ non-woven web should not be used as a surface dressing except for immediate control of bleeding. Avoid packing AVITENE™ tightly in cavities, especially within the bony enclosure of the CNS or within other relatively rigid


Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

cavities where swelling may interfere with normal function or possibly cause necrosis. Avitene™ is not recommended for use in patients sensitive to bovine derived collagen.




Adverse Reactions: The most serious adverse reactions reported which may be related to the use of Avitene™ (MCH) are potentiation of infection including abscess formation, hematoma, wound dehiscence and mediastinitis. Other reported adverse reactions possibly related are adhesion formation, allergic reaction, foreign body reaction and subgaleal seroma (report of a single case). The use of MCH in dental extraction sockets has been reported to increase the incidence of alveolgia. Transient laryngospasm due to aspiration of dry material has been reported following use of MCH in tonsillectomy.

Main technical and functional characteristics

Product Form	Product Code	Characteristics
Avitene™ Microfibrillar Collagen Hemostat Flour	1010010	Mass: 0.5 gram
	1010020	Mass: 1.0 gram
	1010590	Mass: 5.0 grams
Avitene™ Microfibrillar Collagen Hemostat Non-Woven Web	1010080	Size: 35mm x 35mm Basis weight: 270-350 g/m ² Mass: not more than 6 g
	1010090	Size: 70mm x 35mm Basis weight: 270-350 g/m ² Mass: not more than 9 g
	1010110	Size: 70mm x 70mm Basis weight: 270-350 g/m ² Mass: not more than 12 g
Avitene™ Microfibrillar Collagen Hemostat EndoAvitene™	1010150	Mesh size: 50 mm x 15 mm x 1 mm Mass: not more than 3 g Syringe mass: not more than 27 g Syringe diameter: 10 mm Dimpled depth from end of sheath: 0.111 ± 0.040. mm Plunger push force: not more than 26,6893 N (6 pounds) Axial pull force of bonded/welded joints: not less than 44,5 N (10 lb) preloaded sheet diameter: 1 mm usable length: 42 cm Endo-Resulting sheath I.D. after dimpling: 0.275± 0.010 mm
		
		A – 360,93± 6,35 mm (14.21±0.25 inch) B – 20,320± 0,254 mm (0,8±0,01 inch) C – 7,0100+0/- 0,1524 mm(0,276 +0/-0,005 inch)
Sheath		



Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

		
		A – $338,07 \pm 2,54$ mm ($13,31 \pm 0,1$ inch) B – $25,4000 \pm 0,1524$ mm ($1,000 \pm 0,005$ inch) C – $10,0080 \pm 0,0508$ mm ($0,394 \pm 0,002$ inch) D – $8,4836 \pm 0,0508$ mm ($0,334 \pm 0,002$ inch)
	1010260	Mesh size: 50 mm x 5 mm x 1 mm Mass: not more than 0,3 g Syringe diameter: 5 mm Syringe mass: not more than 27 g Dimpled depth from end of sheath: 0.111 ± 0.040 mm Plunger push force: not more than 26,6893 N (6 pounds) Axial pull force of bonded/welded joints: not less than 44,5 N (10 lb) preloaded sheet diameter: 5 mm usable length: 42 cm Endo-Resulting sheath I.D. after dimpling: 0.145
Plunger 		
		A – $20,32 \pm 0,1524$ mm ($0,800 \pm 0,005$ inch) B – $355,6 \pm 0,254$ mm ($14 \pm 0,01$ inch) C – $3,175 \pm 0,1524$ mm ($0,125 \pm 0,005$ inch) D – $3,81 \pm 0,1524$ mm ($0,150 \pm 0,005$ inch) E – $34,036 \pm 0,254$ mm ($1,34 \pm 0,01$ inch) F – $3,2512 \pm 0,1524$ mm ($0,128 \pm 0,005$ inch)
Sheath 		
		A – $338,074 \pm 1,524$ mm ($13,31 \pm 0,06$ inch) B – $299,466 \pm 0,254$ mm ($11,79 \pm 0,01$ inch) C – $38,608 \pm 0,254$ mm ($1,52 \pm 0,01$ inch) D – $25,4000 \pm 0,1524$ mm ($1,000 \pm 0,005$ inch) E – $5,0038 \pm 0,0254$ mm ($0,197 \pm 0,001$ inch) F – $4,2418 \pm 0,0508$ mm ($0,167 \pm 0,002$ inch)
Avitene™ Microfibrillar Collagen Hemostat	1010340	1 gram preloaded flour 2 cm diameter, 16.5 cm usable length

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

SyringeAvitene™	Syringe mass: not more than 55 g
Plunger 	
	A – 203,200 ± 0,254mm (8,00±0,01 inch) B – 18,1356 ± 0,2032 mm (0,714±0,008 inch) C – 15,5956 ± 0,0762 mm (0,614±0,003 inch) D – 17,5260 ± 0,0762 mm (0,690±0,003 inch)
Sheath 	
	A – 203,200 + 0,762/- 0,381 mm (8,00+0,030/-0,015 inch) B – 31,750 ± 0,508 mm (1,25±0,02 inch) C – 1,5748 ± 0,508 mm (0,062±0,02 inch) D – 15,875 ± 0,127 mm (0,625±0,005 inch) E – 3,810 ± 0,127 mm (0,150±0,005 inch) F – 4,826 ± 0,254 mm (0,19±0,01 inch) G – 24,765 ± 0,127 mm (0,975±0,005 inch) H – 19,126 ± 0,127 mm (0,753±0,005 inch)

Tolerance is ± 5% if unless otherwise specified.

Package contents:

Product Form	Product Codes	Units/case
Avitene™ Microfibrillar Collagen Hemostat Flour	1010010	1 unit per package
	1010590	1 unit per package
	1010020	1 unit per package
Avitene™ Microfibrillar Collagen Hemostat Non-Woven Web	1010080	1 unit per package
	1010110	6 unit per package
	1010090	6 unit per package
Avitene™ Microfibrillar Collagen Hemostat EndoAvitene™	1010150	1 unit per package
	1010260	1 unit per package
Avitene™ Microfibrillar Collagen Hemostat SyringeAvitene™	1010340	1 unit per package

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.


Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

Directions for opening pouch:

1. Peel back foil pouch to expose blister pack
2. Using accepted sterile technique, pass blister pack onto sterile pack sterile field
3. Caution: Do not depress syringe plunger prior to positioning of distal end at the surgical site

Dosage And Administration: Avitene™ (MCH) must be applied directly to the source of bleeding. Because of its adhesiveness, it may seal over the exit site of deeper hemorrhage and conceal an underlying hematoma as in penetrating liver wounds. When possible, surfaces to be treated should be compressed with dry sponges immediately prior to application of the dry Avitene™. It is then advantageous to apply pressure over the Avitene™ with a dry sponge for a period of time which varies with the force and severity of bleeding. A minute may suffice for capillary bleeding (e.g., skin graft donor sites, dermatologic curettage) but three to five or more minutes may be required for brisk bleeding (e.g. splenic tears) or high pressure leaks in major artery suture holes. For control of oozing from cancellous bone, it should be firmly packed into the spongy bone surface. After five to ten minutes, excess MCH may be teased away (see Precautions); this can usually be accomplished with blunt forceps and is facilitated by wetting with sterile 0.9% saline solution and irrigation. If breakthrough bleeding occurs in areas of thin application, additional Avitene™ may be applied. The amount required depends, again on the severity of bleeding. In capillary bleeding, one gram will usually be sufficient for a 50 cm² area. Thicker coverage will be required for more brisk bleeding. MCH will adhere to wet gloves, instruments, or tissue surfaces. To facilitate handling, use dry, smooth forceps or a pre-loaded Avitene™ applicator device (SyringeAvitene™ or EndoAvitene™). In neurosurgical and other procedures the non-woven web may conveniently be used by applying small squares to bleeding areas and then covering the sites with moist cottonoid "patties". After five to ten minutes excess MCH may be removed by teasing and irrigation.

AVITENE™ Pre-loaded Applicators

1. Follow accepted sterile technique when removing from foil pouch and introducing into sterile field.
2. Position the distal end of applicator at the bleeding site.
3. Depress the syringe plunger to deliver Avitene™ to the site of application.
4. Excess Avitene™ can be irrigated away.

Usage conditions: The medical device is intended for use at the standard conditions of hospital operating room environment.

Potential consumer: Medical device should only be used by qualified medical personnel with experience in the use of similar products or by medical personnel who have completed training in working with similar products under the guidance of an experienced specialist with experience in the use of similar products. Avitene MCH is a prescription device and as such should be used by a qualified medical professional.

Key Accessories, Medical Devices or Non-Medical Devices intended to be used with Medical Device: There are no accessories used with Avitene™ Microfibrillar Collagen Hemostat.

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

Avitene MCH is not intended for use in conjunction with other medical devices, medicines or other medical technologies. Avitene MCH is applied directly to the source of bleeding by the care provider. The device is not typically positioned with medical imaging. However, Avitene MCH's effect on platelet adhesion and aggregation is not inhibited by heparin.

Pharmaceutical substances: Avitene™ Microfibrillar Collagen Hemostat does not incorporate as an integral part, a substance, which if used separately, may be considered to be a medicinal product as defined in Article 1 of Directive 65/65/EEC and subsequent Directive 2001/83/EC.

Information about materials of animal/human origin: Avitene™ Microfibrillar Collagen Hemostat does not incorporate or utilize any substance or material that is derived from human or animal origin as defined in Article 1.5(f)(g) of Directive 93/42/EEC.

The nature of the starting tissue(s), animal species and geographical source: Avitene is composed of bovine corium. The bovine corium (skin, dermis, hide) raw material is supplied to Davol, Inc. from Devro Pty, Ltd. of Australia. The Bos Taurus cattle are grass fed being born, continuously reared, and slaughtered in Australia.

Davol conducts audits of Devro on a biennial basis to ensure compliance to the documented procedures. Devro certifies compliance of supplied material to ISO 22442-2. The raw material specification specifies certification and traceability requirements.

Animals are slaughtered at the abattoir and transported to Devro for final cleaning, fleshing, cutting, and packaging. All equipment is verified to be clean prior to use. Microbial and pyrogen testing is performed for each lot. Devro closely controls the waste material and removes it from operation efficiently.

Raw hides from the abattoir are transported in covered bins to Devro. After processing at Devro, hides are packed, stored, and transported frozen in 2 food grade bags interleaved with a food grade sheet within a sealed cardboard shipper. Material is maintained frozen for storage and transportation. Each lot is verified upon receipt for acceptability of the required documentation and verification of the chill checkers.

Sterilization methods: Avitene™ Microfibrillar Collagen Hemostat is dry heat sterilized at Davol's Woburn, Massachusetts Operations facility.

Storage and transportation conditions: Medical Device should be stored/ transported at temperatures from -10°C to +35 °C and 65% relative humidity.

Shelf life/Service life: The Avitene™ Microfibrillar Collagen Hemostats are labeled with a 36 month (3 year) shelf life.

Duration of absorption in human body - 8-10 weeks.

Disposal: After use, the applicators are potentially biohazardous products due to exposure to blood and body fluids (class B medical devices under medical wastes classification). Dispose of any potentially biohazardous products in accordance with accepted medical practice and applicable local, state and federal regulations and the institutions written policy for the disposal of potentially biohazardous material.


Unused products (not having contact with blood and/or biological fluids), including those with expired shelf life, belong to class A medical devices under medical waste classification and may be recycled and/or disposed in accordance with applicable national regulations.

Package description: Packaging for the various Avitene™ Microfibrillar Collagen Hemostat product families consists of foil pouches, tyvek lids, blister trays, polycarbonate cup and cartons.

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

Symbols used on packaging:

	Catalogue number
	Sterilized using steam
	CE-mark and identification number of notified body. The product meets the essential requirements of Medical Device Directive 93/42/EEC.
	Batch code
	Single use
	Use-by date
	Manufacturer
	Caution
	Non-pyrogenic
	Non-pyrogenic
	Do not use if package is damaged
	See Instructions for use
	RX only
	Absorbable
	European Authorized Representative for Medical Device Companies

Information about manufacturer:

Legal Manufacturer

Davol, Inc. Subsidiary of C. R. Bard, Inc.

100 Crossings Blvd.

Warwick RI, 02886

USA

Design Center

Davol, Inc. Subsidiary of C. R. Bard, Inc.

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

100 Crossings Blvd.
Warwick RI, 02886
USA

Manufacturing Site

1. Davol Inc. Subsidiary of C. R. Bard, Inc., 100 Crossings Boulevard Warwick, RI 02886, USA
2. Davol, Inc. 160 New Boston Street, Woburn, MA 01801, USA

Warranty: Davol Inc. guarantees the original buyer that this product is free from defects in materials and performance until the expiration date. If it is proved that this product is defective, the buyer can return it to the manufacturer for replacement. Responsibility Davol Inc. subject to the use of this product.

Reclamation: On the territory of the Russian Federation, all questions regarding warranty and product quality may be addressed to the Authorized Representative of the manufacturer.

Authorized Representative in the Russian Federation:

LLC «Bard Rus»
121596, Moscow, Gorbunova st., 2- 204,
+7 499 372-50-02

Avitene, Bard, Davol, EndoAvitene and SyringeAvitene are trademarks and/or registered trademarks of C. R. Bard, Inc. or an affiliate.

Copyright © 1997, 2004, 2009 - 2011 C. R. Bard, Inc. All Rights Reserved.

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