
STATE OF Colorado
COUNTY OF Boulder
SUBSCRIBED AND SWORN TO BEFORE ME
THIS 22 DAY OF June, 2018.
BY Diana Karlin
[Signature]
NOTARY PUBLIC

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Tel: +1 203 845 1000

М.П. / Stamp

"Operating documents" for submission in the Russian Federation:

Jaw Open Sealer/Divider LigaSure, Nano-coated

Manufactured by:

Covidien Llc., USA. 15 Hampshire Street, Mansfield MA 02048, USA

Manufacturing site:

Covidien Llc. USA 5920 Longbow Drive, Boulder, Colorado, 80301, USA

LigaSure Impact™

Curved, Large Jaw, Open Sealer/Divider
Nano-coated

- en** Instructions for Use
- fr** Mode d'emploi
- sv** Bruksanvisning
- da** Brugervejledning
- fi** Käyttöohjeet
- no** Bruksanvisning
- pt** Instruções de uso
- de** Gebrauchsanleitung
- es** Instrucciones de uso
- it** Istruzioni per l'uso
- nl** Gebruiksaanwijzing
- ru** Инструкция по применению
- zh** 使用说明

LF4418**Compatible with:**

FORCETRIAD SW v3.6 or Higher

VLFT10GEN SW v1.1 or Higher

VLLS10GEN SW v1.0 or Higher

LigaSure Impact™

REF LF4418 **Curved, Large Jaw,
Open Sealer/Divider,
Nano-coated**
**36 mm jaw length,
18 cm shaft length**

The LF4418 is designed for use with Covidien electrosurgical generators that include vessel sealing capability. Please refer to the cover page for details on compatible generator models and software versions. If the software version on your generator is lower than required, contact Covidien about software updates.

These instructions assume that the operator is knowledgeable about correct setup and operation of the associated Covidien generator. Refer to the generator user's guide for setup information and for additional warnings and cautions.

The instrument creates a seal by application of radiofrequency (RF) electrosurgical energy to vascular structures (vessels and lymphatics) or tissue bundles interposed between the jaws of the instrument. A blade within the instrument is surgeon actuated to divide tissue.

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Type CF applied part.

Indications for Use

The LigaSure Sealer/Divider is a bipolar electrosurgical instrument intended for use in open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure Sealer/Divider can be used on vessels (arteries, veins, pulmonary arteries, pulmonary veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc.

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This instrument is intended for use ONLY with the Covidien equipment listed on the cover of this document. Use of this instrument with other generators may not result in the desired tissue effect, may result in injury to the patient or surgical team, or may cause damage to the instrument.

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Use the system with caution in the presence of internal or external pacemakers, or other implanted devices. Interference produced by electrosurgical equipment can cause a pacemaker or other device to enter an unsafe mode or permanently damage the device. Consult the device manufacturer or responsible hospital department for further information when use is planned in patients with implanted medical devices.

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
Contact between an active instrument electrode and any metal object (hemostats, staples, clips, retractors, etc.) may increase current flow and may result in unintended surgical effects, such as an effect at an unintended site or insufficient energy deposition.

Precaution

The performance of this single-use device has been tested according to the expected conditions of a single surgical procedure. Subjecting the device to process steps, tools, and chemicals commonly used for reprocessing has been shown to degrade the nano-coating material on the active sealing surfaces. Degradation of the coating can lead to reduced performance, such as increased tissue adhesion.

LF4418 LigaSure Curved, Large Jaw, Open Sealer/Divider

- ① Connector (purple and white)
- ② Rotation wheel (gray)
- ③ Lever (gray and white)
- ④ Activation button (purple)
- ⑤ Cutting trigger (gray)
- ⑥ Handle (gray and white)
- ⑦ Jaws
- ⑧ Shaft
- ⑨ Cable

Precaution

Use caution during surgical cases in which patients exhibit certain types of vascular pathology (atherosclerosis, aneurysmal vessels, etc.). For best results, apply the seal to unaffected vasculature.

Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages.

Setup

Warning

Electric Shock Hazard Do not connect wet accessories to the generator.

Confirm proper energy platform settings before proceeding with surgery.

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Position instrument cords to avoid contact with the patient or other cords. Do not wrap cords around metal objects. This may induce currents that can lead to shocks, fires, or injury to the patient or surgical team.

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Do not use in the presence of flammable anesthetics or oxidizing gases (such as nitrous oxide (N₂O) and oxygen) or in close proximity to volatile solvents (such as ether or alcohol) as explosion may occur.

Because of concerns about the carcinogenic and infectious potential of electrosurgical by-products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke evacuation equipment should be used.

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1. Remove instrument from tray by firmly pulling on the handle (6). Do not pull on the instrument's jaws (7) or cable (9).
2. Insert the connector (1) into the receptacle on the generator. Follow the instructions in the generator user's guide to complete the setup procedure.

During Surgery

Warning

Fire Hazard Do not place instruments near or in contact with flammable materials (such as gauze or surgical drapes). Instruments that are activated or hot from use may cause a fire. When not using instruments, place them in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.

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Place the vessel or tissue in the center of the jaws. To avoid incomplete sealing, do not grasp tissue beyond the electrode surface; do not place tissue in the jaw hinge.

Do not use this instrument on vessels in excess of 7 mm in diameter.

Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to the instrument may carry electrical current or heat. This may cause unintended burns to the patient. Aspirate fluid from around the instrument jaws before activating the instrument.

Always keep the external surface of the instrument jaws away from adjacent tissue while activating the instrument.

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Do not activate the LigaSure system until the lever has been properly latched. Activating the device before this is done may result in improper sealing and may increase thermal spread to tissue outside the surgical site.

Do not activate the instrument while in contact with or near other instruments, as localized burns to the patient or physician may occur.

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Eliminate tension on the tissue while sealing and cutting to ensure proper function.

Do not attempt to seal or cut over clips or staples as incomplete seals/damage to the cutting blade will occur. Contact between an active electrode and any metal objects may result in alternate site burns or incomplete seals.

Use caution during surgical cases in which patients exhibit incompressible tissue. There may be limitations associated with effective use of the device in these situations.

Although the jaws of the instrument will accommodate tissue greater than 15 mm thick, exceeding the tissue thickness limit of 15 mm could result in compromised seals.

Precaution

Activate the instrument *only* when ready to deliver electrosurgical current and when the active tip is in view.

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Rotating the Instrument Jaws

Notice

Do not turn the rotation wheel (2) when the lever (3) is latched. Product damage may occur.

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Precaution

Keep the instrument jaws clean. Build-up of eschar may reduce sealing and/or cutting effectiveness. Wipe jaw surfaces and edges with a sterile, wet gauze pad as needed.

1. Open the jaws by pushing forward on the lever.
2. Grasp the intended vessel and/or tissue in the center of the jaws.
3. Squeeze the lever until it latches in place.
4. Relax the hand on the lever once latched. The device will maintain the latched position.

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5. Activate the instrument by either of the following methods:
 - Press and hold the activation button (4) on the back of the instrument.
 - By using an appropriate footswitch associated with the generator.

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Note: If both the activation button and footswitch pedal are activated during the same seal cycle, the generator delivers energy from the activation source it detects first.

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
6. Release the activation button on the instrument, or the footswitch pedal when the seal cycle is complete.
7. To seal adjacent tissue, overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin. Failure to overlap seals when a second seal is desired may compromise seal integrity.

Cutting Tissue

Warning

Do not overfill the jaws of the instrument with tissue. This may damage the cutting mechanism or cause the blade to deploy outside of its guiding features, possibly resulting in difficulty opening the jaws or unintended injury to the user or patient.

Energy-based devices, such as ESU pencils or ultrasonic scalpels, that are associated with thermal spread should not be used to transect seals.

① < 2 mm – Confirm that the jaws are in the closed position prior to activating the cutter.

Notice

Do not engage the cutting mechanism over sutures, clips, staples, or other metal objects as damage to the cutter may occur.

1. To engage the cutting mechanism, pull the cutting trigger (5) completely back toward the body of the instrument. The instrument blade transects the seal an average of 2 mm from the distal end of the seal plate.

- Open the jaws by squeezing the lever until it unlocks, then push the lever completely forward.

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Cleaning the Instrument During Use

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Do not activate the instrument or the cutting trigger while cleaning the jaws. Injury to operating room personnel may result.

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Troubleshooting information	
The following is a list of troubleshooting suggestions for situations encountered when using the instrument with compatible Covidien vessel sealing generators. For details on specific situations, refer to the corresponding generator user's guide or the generator quick reference guide.	
Alert situations	When an alert condition occurs, energy delivery stops, the generator produces a sequence of pulsed tones, and an alert will be displayed on the generator. Do Not Cut the Vessel. The user should inspect the seal site and instrument before proceeding. After the alert condition has been corrected, energy delivery will be immediately available.
Troubleshooting steps	1) Release the footswitch pedal or activation button, if still engaged. 2) Open the instrument jaws and inspect for a successful seal. 3) Follow the suggested corrective actions on the generator screen, the generator quick reference card, or in the generator user's guide. 4) If possible, reposition the instrument and regrasp tissue in a location that overlaps the previous seal, then reactivate the seal cycle.
Reasons for alert	Too little tissue between the jaws – The user is grasping thin tissue or not enough tissue; open the jaws and confirm that a sufficient amount of tissue is inside the jaws. If necessary, increase the thickness of tissue that is grasped and reactivate the seal cycle. Too much tissue between the jaws – The user is grasping too much tissue; open the jaws, reduce the amount of tissue that is grasped, and reactivate the seal cycle. Activating on a metal object – Avoid grasping objects, such as staples, clips, or encapsulated sutures in the jaws of the instrument. Dirty jaws – Use a wet gauze pad to clean surfaces and edges of instrument jaws. Excess Fluids in the Surgical Field – Minimize or remove excess fluids from around the instrument jaws. Activation switch released before seal complete tone – The footswitch or activation button was released before the seal cycle was complete. Maximum seal cycle time has been reached – The system needs more time and energy to complete the seal cycle.

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Product performance of the device was established in a chronic in-vivo porcine model. The results showed that no animals studied experienced any hemostatic complications related to the device during the 21-day survival period. A variety of tissue types and vessels were evaluated to demonstrate effective sealing in arteries and veins up to and including 7 mm.

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Vessel type	Vessel identification	Vessel size range
A/V Bundle	Splenic Mesentery	< 1.5 mm bundles
	Gastrosplenic	3.5 mm – 4.5 mm arteries within bundles
	Short Gastric	4.5 mm – 5.0 mm bundles
	Ovarian Pedicle	Bundles up to 4.0 mm
	Broad Ligament	Bundles up to 4.0 mm
Artery	Splenic	5.0 mm – 7.0 mm
	Renal	4.0 mm – 7.0 mm
Vein	Splenic	7.0 mm
	Renal	4.5 mm – 6.5 mm

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For information regarding California Proposition 65, please refer to www.covidien.com/caprop65.

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O: indicates that the content of the toxic and hazardous substance in all the Homogenous Materials of the part is below the concentration limit requirement as described in GB/T 26572.

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REV 11/2016

LigaSure™

Blunt Tip Sealer/Divider
Nano-coated

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REF LF1837 **Blunt Tip Laparoscopic Sealer/Divider**
Nano-coated
5 mm – 37 cm

REF LF1844 **Blunt Tip Laparoscopic Sealer/Divider**
Nano-coated
5 mm – 44 cm

The LigaSure Blunt Tip Sealer/Divider is designed for use with Covidien electrosurgical generators that include vessel sealing capability. Please refer to the cover page for details on compatible generator models and software versions. If the software version on your generator is lower than required, contact Covidien about software upgrades.

These instructions assume that the operator is knowledgeable about correct set-up and operation of the associated Covidien generator. Refer to the generator user's guide for set up information and for additional warnings and precautions.

The instrument creates a seal by application of radiofrequency (RF) electrosurgical energy to vascular structures (vessels and lymphatics) or tissue bundles interposed between the jaws of the instrument. A blade within the instrument is surgeon actuated to divide tissue. Shorter shaft lengths are typically used for open procedures while longer shaft lengths are typically used for minimally invasive procedures.

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Type CF applied part

Indications for Use

The LigaSure Sealer/Divider is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue

bundles, and lymphatics is desired. The LigaSure Sealer/Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. Procedures may include, but are not limited to, Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc.

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
The instrument is intended for use ONLY with the Covidien equipment listed on the cover of this document. Use of this instrument with other generators may not result in the desired tissue effect, may result in injury to the patient or surgical team, or may cause damage to the instrument.

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When this instrument is used with an energized endoscope, the leakage current from the instrument and the endoscope are additive. The patient may be exposed to unexpected levels of leakage current if this instrument is used with an energized endoscope that is not type CF.

LF1823, LF1837, and LF1844 LigaSure Blunt Tip Sealer/Divider

- ① Jaws
- ② Shaft
- ③ Rotation wheel (gray)
- ④ Handle (gray and white)
- ⑤ Activation button (purple)
- ⑥ Connector (purple and white)
- ⑦ Cable
- ⑧ Lever (gray and white)
- ⑨ Cutting trigger (gray)

Precaution

The performance of this single-use device has been tested according to the expected **conditions of a single surgical procedure**. Subjecting the device to process steps, tools, and chemicals commonly used for reprocessing has been shown to degrade the nano-coating material on the active sealing surfaces. Degradation of the coating can lead to reduced performance, such as increased tissue adhesion.

Warning

In minimally invasive surgery, inspect the outer surfaces of the instrument before insertion through the cannula to ensure that there are no rough or sharp edges that could damage tissue.

Contact between an active instrument electrode and any metal objects (hemostats, staples, clips, retractors, etc.) may increase current flow and may result in unintended surgical effects, such as an effect at an unintended site or insufficient energy deposition.

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Position instrument cords to avoid contact with the patient or other cords. Do not wrap cords around metal objects. This may induce currents that could lead to shocks, fires, or injury to the patient or surgical team.

Examine all LigaSure system and instrument connections before using. Improper connections may result in arcing, sparks, accessory malfunction, or unintended surgical effects.

Inspect the instrument and cords for breaks, cracks, nicks, or other damage before use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team or cause damage to the instrument. If damaged, do not use.

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Do not use in the presence of flammable anesthetics or oxidizing gases, such as nitrous oxide (N₂O) and oxygen, or in close proximity to volatile solvents (such as ether or alcohol) as explosion may occur.

Because of concerns about the carcinogenic and infectious potential of electrosurgical by-products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke evacuation equipment should be used in both open and minimally invasive procedures.

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1. Remove instrument from tray by firmly pulling on the handle (4). Do not pull on the instrument's jaws (1) or cable (7).
2. Insert the connector (6) into the receptacle on the generator. Follow the instructions in the generator user's guide to complete the setup procedure.

During Surgery**Tissue Manipulation and Dissection**

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Keep the cord free from the jaw and latch area of the instrument.

Warning

Fire Hazard– Do not place instruments near or in contact with flammable materials (such as gauze, surgical drapes, or flammable gases). Instruments that are activated or hot from use may cause a fire. When not using instruments, place them in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.

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Rotating the Instrument Jaws**Notice**

Do not turn the rotation wheel (3) when the lever (8) is latched. Product damage may occur.

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1. Open the jaws by pushing forward on the lever.
2. Grasp the intended vessel and/or tissue in the center of the jaws.

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Use caution when grasping, manipulating, sealing, and dividing large tissue bundles.

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Inadvertent activation or movement of the activated instrument outside of the field of vision may result in injury to the patient.

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3. Squeeze the lever until it latches in place.

Precaution

Do not activate the LigaSure system until the lever has been latched. Activating the system before latching may result in improper sealing and may increase thermal spread to tissue outside of the surgical site.

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Keep the instrument jaws (1) clean. Build-up of eschar may reduce the seal and/or cutting effectiveness. Wipe jaw surfaces and edges with a wet gauze pad as needed.

4. Relax the hand on the lever once latched. The device will maintain the latched position.

Notice

Do not apply additional hand force to the lever during sealing to ensure proper function.

5. The LigaSure system can be activated by either of the following two methods:
- Press and hold the activation button (5) on the back of the instrument.
 - or
 - By using an appropriate footswitch associated with the generator.

Refer to the generator user's guide for footswitch instructions.

Note: If both the activation button and footswitch pedal are activated during the same seal cycle, the vessel sealing system delivers energy from the activation source it detects first.

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A tone with multiple pulses indicates that the seal cycle was not completed. Refer to the Troubleshooting section on page 6 for possible causes and corrective actions. Do not cut tissue until you have verified that there is an adequate seal.

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6. Release the activation button on the instrument or the footswitch pedal when the seal cycle is complete.
7. To seal adjacent tissue, overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.

Cutting Tissue

Warning

Prior to cutting the seal, inspect the vessel or tissue to ensure proper sealing.

Energy-based devices, such as ES pencils or ultrasonic scalpels that are associated with thermal spread, should not be used to transect seals.

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1. To engage the cutting mechanism, pull the cutting trigger (9) completely back toward the body of the instrument.

Important

If the cutting trigger does not automatically return to position, unlatch and open the lever to manually return the cutting trigger.

2. Open the jaws by squeezing the lever until it unlocks, then push it completely forward.

Warning

If the lever cannot be unlatched following use, open the device by forcing the lever forward from the handle. The device will no longer function properly and must be discarded.

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Troubleshooting Information	
The following is a list of troubleshooting suggestions for situations encountered when using the instrument with compatible Covidien vessel sealing generators. For details on specific situations, refer to the corresponding generator user's guide or the generator quick reference guide.	
Alert Situations	When an alert condition occurs, energy delivery stops, the generator produces a sequence of pulsed tones, and an alert will be displayed on the generator. Do Not Cut the Vessel. The user should inspect the seal site and instrument before proceeding. After the alert condition has been corrected, energy delivery will be immediately available.
Troubleshooting Steps	1) Release the footswitch pedal or activation button, if still engaged. 2) Open the instrument jaws and inspect for a successful seal. 3) Follow the suggested corrective actions on the generator screen, the generator quick reference card, or in the generator user's guide. 4) If possible, reposition the instrument and regrip tissue in a location that overlaps the previous seal, then reactivate the seal cycle.

Troubleshooting Information	
Reasons for Alert	Too little tissue between the jaws – The user is grasping thin tissue or not enough tissue; open the jaws and confirm that a sufficient amount of tissue is inside the jaws. If necessary, increase the thickness of tissue that is grasped and reactivate the seal cycle. Too much tissue between the jaws – The user is grasping too much tissue; open the jaws, reduce the amount of tissue that is grasped, and reactivate the seal cycle. Activating on a metal object – Avoid grasping objects, such as staples, clips, or encapsulated sutures in the jaws of the instrument. Dirty jaws – Use a wet gauze pad to clean surfaces and edges of instrument jaws. Excess Fluids in the Surgical Field – Minimize or remove excess fluids from around the instrument jaws. Activation switch released before seal complete tone – The footswitch or activation button was released before the seal cycle was complete. Maximum seal cycle time has been reached – The system needs more time and energy to complete the seal cycle.

After Surgery

Discard the instrument after use according to hospital policy for biohazards and sharps. **Do not resterilize.**

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Clearance of the device in the United States was not based on human clinical testing; instead, performance of the device was established in a chronic in vivo porcine model. The results showed that no animals studied experienced any hemostatic complications related to the device during the 21-day survival period. The evaluation was performed on a variety of tissue types and vessels, up to and including 7 mm, to demonstrate effective sealing.

Vessel Type	Tissue/Vessel Name	Size Range
A/V Bundle	Omentum	<1.0 mm bundles
	Gastrosplenic	4.0 – 4.5 mm arteries within bundles
	Ovarian Pedicle	Bundles up to 5.0 mm
	Broad Ligament	Bundles up to 4.0 mm
	Short Gastric	4.5– 7.0 mm bundles
Isolated Vessel (Artery, Vein)	Renal	4.5 mm – 7.0 mm
	Splenic	5.5 mm – 7.0 mm

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California Proposition 65 Statement

For information regarding California Proposition 65, please refer to www.covidiem.com/caprop65.

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X: indicates that the content of the toxic and hazardous substance in at least one Homogenous Material of the part exceeds the concentration limit requirement as described in GB/T 26572 (Information provided in this table is solely for compliance with SJ/T 11364 and the Management Measures for the Control of Pollution by Electronic Information Products.)

STERILE EO

Rx
ONLY

CE
0086

Part No. PT00016000

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM brands are trademarks of their respective owner.

© 2016 Covidien.

 Covidien Ilc,
15 Hampshire Street, Mansfield, MA 02048 USA.
 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Made in USA. Printed in USA.

www.covidien.com

REV 10/2016

LigaSure™

Maryland Jaw Sealer/Divider
One-step Sealing, Nano-coated

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VLLS10GEN SW v1.0 or Higher

LigaSure™

- [REF] LF1923 Maryland Jaw Open Sealer/Divider, One-step Sealing Nano-coated 5 mm – 23 cm**
- [REF] LF1937 Maryland Jaw Laparoscopic Sealer/Divider, One-step Sealing Nano-coated 5 mm – 37 cm**
- [REF] LF1944 Maryland Jaw Laparoscopic Sealer/Divider, One-step Sealing Nano-coated 5 mm – 44 cm**

The Maryland Jaw Sealer/Divider is designed for use with Covidien electrosurgical generators that include vessel sealing capability. Please refer to the cover page for details on compatible generator models and software versions. If the software version on your generator is lower than required, contact Covidien about software updates at <http://www.medtronic.com/covidien/products/electrosurgical-hardware>.

These instructions assume that the operator is knowledgeable about correct set-up and operation of the associated Covidien generator. Refer to the generator user's guide for set-up information and for additional warnings and precautions.

The instrument creates a seal by application of radiofrequency (RF) electrosurgical energy to vascular structures (vessels and lymphatics) or tissue bundles interposed between the jaws of the instrument. A blade within the instrument is surgeon-actuated to divide tissue. Shorter shaft lengths are typically used for open procedures while longer shaft lengths are typically used for laparoscopic procedures.

Maximum rated voltage: 288 V_{peak}

Not made with natural rubber latex

Do not use if package is opened or damaged

Type CF applied part

Indications for Use

The LigaSure Sealer/Divider is a bipolar electrosurgical instrument intended for use in minimally invasive or open surgical procedures where ligation and division of vessels, tissue bundles, and lymphatics is desired. The LigaSure Sealer/Divider can be used on vessels (arteries and veins) up to and including 7 mm. It is indicated for use in general surgery and such surgical specialties as urologic, vascular, thoracic, and gynecologic. These may include, but are not limited to, such procedures as Nissen fundoplication, colectomy, cholecystectomy, adhesiolysis, hysterectomy, oophorectomy, etc.

The LigaSure system has not been shown to be effective for tubal sterilization or tubal coagulation for sterilization procedures. Do not use the LigaSure system for these procedures.

General Warnings**Warning**

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices without appropriate regulatory authorization may result in bio-incompatibility, infection, or product failure risks to the patient.

These instruments are intended for use ONLY with the Covidien equipment listed on the cover of this document. Use of this instrument with other generators may not result in the desired tissue effect, may result in injury to the patient or surgical team, or may cause damage to the instrument.

Do not use the LigaSure system unless properly trained to use it in the specific procedure being undertaken. Use of this equipment without such training may result in serious unintended patient injury.

Use the system with caution in the presence of internal or external pacemakers or other implanted devices. Interference produced by electrosurgical equipment can cause a pacemaker or other device to enter an unsafe mode or permanently damage the device. Consult the device manufacturer or responsible hospital department for further information when use is planned in patients with implanted medical devices.

When this instrument is used with an energized endoscope, the leakage current from the instrument and the endoscope are additive. The patient may be exposed to unexpected levels of leakage current if this instrument is used with an energized endoscope that is not a type CF applied part.

Warning

In minimally invasive surgery, inspect the outer surfaces of the instrument before insertion through the cannula to ensure there are no rough or sharp edges to damage tissue.


Contact between an active instrument electrode and any metal objects (hemostats, staples, clips, retractors, etc.) may increase current flow and may result in unintended surgical effects or insufficient energy deposition.

The safe and effective use of RF energy depends on many factors solely under the control of the operator. There is no substitute for properly trained and vigilant personnel. It is important that the operating instructions supplied with this or any other medical equipment be read, understood, and followed.

Precaution

Use caution during surgical procedures in which patients exhibit certain types of vascular pathology (atherosclerosis, aneurysmal vessels, etc.). For best results, apply the seal to unaffected vasculature.

The performance of this single-use device has been tested according to the expected conditions of a single surgical procedure. Subjecting the device to process steps, tools, and/or chemicals commonly used by third-party re-processors may negatively affect its performance. These have the potential to degrade the LigaSure nano-coating technology on the sealing surfaces, which may lead to increased tissue adhesion.

LF1923, LF1937, LF1944 LigaSure Maryland Jaw Sealer/Divider

- ① Jaws
- ② Shaft
- ③ Rotation wheel (gray)
- ④ Handle (gray and white)
- ⑤ Activation button (purple)

- ⑥ Connector
- ⑦ Cable
- ⑧ Lever (gray and white)
- ⑨ Cutting trigger (gray)

Setup

Warning

Electric Shock Hazard – Do not connect wet accessories to the LigaSure system.

Position instrument cords to avoid contact with the patient or other cords. Do not wrap cords around metal objects. This may induce currents that could lead to shock, fire, or injury to the patient or surgical team.

Warning

Examine all LigaSure system and instrument connections before using. Improper connection may result in arcing, sparks, accessory malfunction, or unintended surgical effects.

Inspect the instrument and cords for breaks, cracks, nicks, or other damage before use. Failure to observe this caution may result in injury or electrical shock to the patient or surgical team, or cause damage to the instrument. If damaged, do not use.

Confirm proper LigaSure system settings before proceeding with surgery.

Do not use in the presence of flammable anesthetics or oxidizing gases, such as nitrous oxide (N₂O) and oxygen, or in close proximity to volatile solvents (such as ether or alcohol) as an explosion may occur.

Because of concerns about the carcinogenic and infectious potential of electrosurgical by-products (such as tissue smoke plume and aerosols), protective eye wear, filtration masks, and effective smoke-evacuation equipment should be used in both open and minimally invasive procedures.

Connect adaptors and accessories to the electrosurgical unit only when the unit is off or in standby mode. Failure to do so may result in injury or electrical shock to the patient or operating personnel.

Precaution

Inspect packaging for damage. If damaged, do not use.

If the generator provides multiple power settings, use the lowest power needed to achieve the intended effect.

1. Remove instrument from tray by firmly pulling on the handle (4). Do not pull on the instrument's jaws (1) or cable (7).
2. Insert the connector (6) into the receptacle on the generator. Follow the instructions in

the generator user's guide to complete the setup procedure.

During Surgery

Tissue Manipulation and Dissection

The instrument can be used to manipulate and dissect tissue with the jaws either open or closed.

Warning

Avoid placing fingers between the lever and the handle, or between the lever and the trigger, or in the jaws. Injury to the user may result.

Use caution when handling the instrument between uses to avoid accidental activation of the LigaSure system. Do not place the instrument on the patient or drapes when not in use.

Keep the cord free from the jaw and lever of the instrument.

Fire Hazard – Do not place instruments near or in contact with flammable materials (such as gauze, surgical drapes, or flammable gases). Instruments that are activated or hot from use may cause a fire. When not using instruments, place them in a clean, dry, highly visible area not in contact with the patient. Inadvertent contact with the patient may result in burns.

For minimally invasive procedures, be alert to these potential hazards:

- Do not use hybrid trocars that are comprised of both metal and plastic components. Capacitive coupling of RF current may cause unintended burns.
- Use the appropriately sized trocar to allow for easy insertion and extraction of the instrument.
- Carefully insert and withdraw the instrument through the cannula to avoid damage to the device and/or injury to the patient.
- Close jaws using device lever before insertion/extraction in the trocar.

Rotating the Instrument Jaws

Notice

Do not turn the rotation wheel (3) when the lever (8) is fully depressed. Product damage may occur.

Turn the rotation wheel on the instrument until the jaws are in the required position.

Grasping and Manipulating Tissue

To grasp tissue with the device, place the tissue in the jaws and pull back on the handle. The first click indicates the end of the grasp zone and alerts the user that additional pressure will activate energy.

Precaution

When grasping or manipulating tissue without intending to activate the device, avoid excess pressure to the lever. Additional pressure after the first click will depress the activation button (5) and deliver energy.

Sealing Vessels and Tissue Bundles

Warning

Do not use this instrument on vessels larger than 7 mm in diameter.

If the instrument shaft is visibly bent, discard and replace the instrument. A bent shaft may prevent the instrument from sealing or cutting properly.

Do not place the vessel and/or tissue in the jaw hinge. Place the vessel and/or tissue in the center of the jaws.

Conductive fluids (e.g., blood or saline) in direct contact with or in close proximity to the instrument may carry electrical current or heat, which may cause unintended burns to the patient. Aspirate fluid from around the instrument jaws before activating the instrument.

Keep the external surface of the instrument jaws away from adjacent tissue while activating the LigaSure system or unintended injury may result.

During a seal cycle, energy is applied to the tissue between the instrument jaws. This energy may cause water to be converted into steam. The thermal energy of steam may cause unintended injury to tissues in close proximity to the jaws. Care should be taken in surgical procedures occurring in confined spaces in anticipation of this possibility.

Eliminate tension on the tissue when sealing and cutting to ensure proper function.

Use caution when grasping, manipulating, sealing, and dividing tissue bundles.

Do not bend instrument shaft.

Warning

Do not attempt to seal or cut over clips or staples as incomplete seals will be formed. Contact between an active electrode and any metal objects may result in alternate site burns or incomplete seals.

The surfaces of the jaws may remain hot enough to cause burns after the RF current is deactivated.

Inadvertent activation or movement of the activated instrument outside of the field of vision may result in injury to the patient or surgical team.

Do not activate the LigaSure system in an open-circuit condition. Activate the system only when the instrument is in direct contact with the target tissue to lessen the possibility of unintended burns.

Do not activate the instrument while instrument jaws are in contact with, or in close proximity to, other instruments including metal cannulas, as localized burns to the patient or physician may occur.

Precaution

Do not overfill the jaws of the instrument with tissue, as this may reduce device performance.

Keep the instrument jaws clean. Build-up of eschar may reduce the seal and/or cutting effectiveness. Wipe jaw surfaces and edges with a wet gauze pad as needed.

Sealing with Hand Activation

Notice

Closing the lever to the second click activates or reactivates energy delivery if hand-activation is being used.

1. Ensure hand-activation is enabled on the appropriate receptacle. Refer to the generator user's guide if needed.
2. Open the jaws by pushing forward on the lever.
3. Squeeze the lever to grasp the intended vessel and/or tissue in the center of the jaws. The first click indicates the user has reached the end of the grasp zone.
4. Squeeze the lever until the button clicks a second time. Continue holding the lever closed. At the second click, energy is delivered. Continue to hold the lever closed until the seal cycle is complete.

Precaution

The lever must be continually held with the activation button fully depressed until the seal cycle is complete. The lever does not latch into the activation position.

A continuous tone sounds to indicate the activation of RF energy. When the activation cycle is complete, a two-pulsed Seal-Cycle-Complete tone sounds and RF output ceases.

Notice

The surgeon may inspect the seal before cutting the vessel or tissue. After inspecting the seal, the surgeon should create a second seal adjacent to the first seal before cutting, as described below.

A tone with multiple pulses indicates that the seal cycle was not completed. Refer to the Troubleshooting section on page 7 for possible causes and corrective actions. Do not cut tissue until you have verified that there is an adequate seal.

5. Open the jaws to release tissue.
6. To seal adjacent tissue, overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.

Notice

Keep the activation button dry and clean.

Sealing with Footswitch Activation

A footswitch can be used instead of the activation button. Ensure that the footswitch is connected to the footswitch receptacle that corresponds to the instrument in use.

Warning

Activating energy delivery with a footswitch when the activation button is not fully depressed may result in improper sealing and increase thermal spread to tissue outside the surgical site. Proper pressure is being applied to the tissue when the lever keeps the activation button fully depressed.

1. Ensure hand-activation is disabled on the appropriate port.
2. Squeeze the lever to grasp the intended vessel and/or tissue in the center of the jaws. The first click will indicate the end of the grasp zone has been reached.
3. Squeeze the lever until the button clicks a second time. Continue holding the lever closed until the seal cycle is complete.

4. Press and hold the footswitch pedal to activate energy until the seal cycle is complete.

A continuous tone sounds to indicate the activation of RF energy. When the activation cycle is complete, a two-pulsed Seal-Cycle-Complete tone sounds and RF output ceases.

Notice

The surgeon may inspect the seal before cutting the vessel or tissue. After inspecting the seal, the surgeon should create a second seal adjacent to the first seal before cutting, as described below.

A tone with multiple pulses indicates that the seal cycle was not completed. Refer to the Troubleshooting section on page 7 for possible causes and corrective actions. Do not cut tissue until you have verified that there is an adequate seal.

5. To seal adjacent tissue, overlap the edge of the existing seal. The second seal should be distal to the first seal to increase seal margin.

Notice

Keep the activation button dry and clean.

Cutting Tissue

Warning

Energy-based devices, such as electrosurgical pencils or ultrasonic scalpels that are associated with thermal spread, should not be used to transect seals.

Notice

Do not engage the cutting mechanism over clips, staples, or other metal objects as damage to the cutter may occur.

1. To activate the cutting mechanism:
 - Grasp the tissue firmly in the jaws by applying steady pressure on the lever.
 - Pull the cutting trigger (9).
 - Release the cutting trigger to retract the cutting blade.

Precaution

Failure to maintain steady pressure on the lever while cutting can result in inadvertent reactivation of energy.

Precaution

If the cutting trigger does not automatically return to position, open the lever to manually return the cutting trigger.

- 2. Open the jaws by pushing forward on the lever to release tissue.

Notice

Keep the activation button dry and clean.

Cleaning the Instrument during Use

Warning

Inspect the instrument jaws prior to cleaning to ensure the blade is not deployed.

Do not activate the instrument or cutting trigger while cleaning the jaws. Injury to operating room personnel may result.

Wipe jaw surfaces and edges with a wet gauze pad as needed.

Notice

Do not attempt to clean the instrument jaws by activating the instrument on wet gauze. Product damage may occur.

Remove any embedded tissue from blade track and jaw hinge area.

Do not clean the instrument jaws with a scratch pad or other abrasives.

Troubleshooting

The following is a list of troubleshooting suggestions for situations encountered when using the instrument with compatible Covidien vessel sealing generators. For details on specific situations, refer to the corresponding generator user's guide or the generator quick reference guide.

Alert Situations

When an alert condition occurs, energy delivery stops. After the alert condition has been corrected, energy delivery will be immediately available.

Troubleshooting information	
The following is a list of troubleshooting suggestions for situations encountered when using the instrument with compatible Covidien vessel sealing generators. For details on specific situations, refer to the corresponding generator user's guide or the generator quick reference guide.	
Alert situations	When an alert condition occurs, energy delivery stops, the generator produces a sequence of pulsed tones, and an alert will be displayed on the generator. Do Not Cut the Vessel. The user should inspect the seal site and instrument before proceeding. After the alert condition has been corrected, energy delivery will be immediately available.
Troubleshooting steps	1) Release the footswitch pedal or activation button, if still engaged. 2) Open the instrument jaws and inspect for a successful seal. 3) Follow the suggested corrective actions on the generator screen, the generator quick reference card, or in the generator user's guide. 4) If possible, reposition the instrument and regrasp tissue in a location that overlaps the previous seal, then reactivate the seal cycle.

	Troubleshooting information
Reasons for alert	Too little tissue between the jaws – The user is grasping thin tissue or not enough tissue; open the jaws and confirm that a sufficient amount of tissue is inside the jaws. If necessary, increase the thickness of tissue that is grasped and reactivate the seal cycle. Too much tissue between the jaws – The user is grasping too much tissue; open the jaws, reduce the amount of tissue that is grasped, and reactivate the seal cycle. Activating on a metal object – Avoid grasping objects, such as staples, clips, or encapsulated sutures in the jaws of the instrument. Dirty jaws – Use a wet gauze pad to clean surfaces and edges of instrument jaws. Excess Fluids in the Surgical Field – Minimize or remove excess fluids from around the instrument jaws. Activation switch released before seal complete tone – The footswitch or activation button was released before the seal cycle was complete. Maximum seal cycle time has been reached – The system needs more time and energy to complete the seal cycle.

After Surgery

Discard the instrument after use according to the facility's policy for biohazards and sharps.
Do not resterilize.

Pre-Clinical Study

Notice

There is no animal data qualified to predict the effectiveness of this device in sealing vessels containing atherosclerotic plaque.

Product performance of the device was established in a chronic in-vivo porcine model. The results showed that no animals studied experienced any hemostatic complications related to the device during the 21-day survival period. A variety of tissue types and vessels was evaluated to demonstrate effective sealing in arteries and veins up to and including 7 mm.

The United States clearance of this device was not based on human clinical testing.

Vessel Type	Tissue/Vessel Name	Vessel Size Range
A/V bundle	Splenic Mesentery	≤ 2.0 mm
	Gastrosplenic	3.0 mm – 4.5 mm arteries within bundles
	Ovarian Pedicle	Bundles up to 5.0 mm
	Broad Ligament	Bundles up to 4.0 mm
	Short Gastric	4.0 mm – 6.0 mm bundles
Artery	Renal	3.5 mm – 7.0 mm
	Splenic	4.5 mm – 7.0 mm
Vein	Renal	3.0 mm – 7.0 mm
	Splenic	7.0 mm

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O: Indicates that the content of the toxic and hazardous substance in all the Homogenous Materials of the part is below the concentration limit requirement as described in GB/T 26572.

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Part No. PT00016037

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM* brands are trademarks of their respective owner.

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