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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of XU RONGCHUN of YANGZHOU FARTLEY MEDICAL INSTRUMENT TECHNOLOGY CO., LTD. and the seal of the said company on the annexed DOCUMENT are genuine.

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书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

Yangzhou Fartley Medical Instrument
Technology Co., Ltd
Beizhou Road, Lidian Town, Guangling District,
Yangzhou, P.R.C. 225106

Tel: 0514-87261752
0514-89888084
Fax: 0514-87268752
Wed: www.fartley.com
E-mail: sales@fartley.com

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To be submitted to the Federal Service for Surveillance in
Healthcare (Roszdravnadzor) in Russian Federation

109074, Moscow, Slavyanskaya Square, 4/1

Date: 26.12.2023

We, Yangzhou Fartley Medical Instrument
Technology Co., Ltd., located at Beizhou Road, Lidian
Town, Guangling District, Yangzhou 225106, Jiangsu P.R.
China, established and reputable manufacturer of the
medical device "Disposable Endoscope Injection
Needle", hereby state that the attached documents are true:

1. Instruction manual for medical device
«Disposable Endoscope Injection Needle».


2 Addendum to Instruction manual for medical
device «Disposable Endoscope Injection Needle».

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
Disposable Endoscope Injection Needles

ECC^{INJECT}

INSTRUCTION MANUAL

EN	Disposable Endoscope Injection Needles	3
DE	Endoskopische Einweg-Injektionsnadeln	5
FR	Aiguilles à injection jetables pour endoscope	7
NL	Endoscopische injectienaalden voor eenmalig gebruik	9
NO	Engangs injeksjonsnåler for endoskop	11
SE	Kasserbara injektionsnålar för endoskop	13
FI	Kertakäyttöiset endoskooppi-injektioneulat	15
DK	Endoskopkanyler til engangsbrug	17
EE	Ühekordselt kasutatavad endoskoobi süstimisnõelad	19
LV	Vienreizējās lietošanas endoskopa injekcijas adatas	21
LT	Vienkartinės endoskopinės injekcinės adatos	23
RU	Одноразовые эндоскопические инъекционные иглы	25
KZ	Бір рет пайдаланылатын эндоскопты инъекциялау ижелері	27
CZ	Jednorázové injekční jehly pro endoskop	29
PL	Jednorazowe igły iniekcyjne do stosowania z endoskopem	31
SK	Jednorázové endoskopické injekčné ihly	33
UA	Ін'єкційні ендоскопічні голки одноразового використання	35
RO	Ace de injectare endoscopică consumabile	37
HU	Egyszer használatos injekciós tűk endoszkóphoz	39
BG	Еднократни ендоскопски инжекционни игли	41
SI	Endoskopske injekcijske igle za enkratno uporabo	43
HR	Jednokratne igle za endoskopsko ubrizgavanje	45
SR	Igle za injekciju za endoskop za jednokratnu upotrebu	47
IT	Aghi monouso per iniezione endoscopica	49
ES	Agujas de inyección de endoscopio desechables	51
PT	Agulhas de injeção para endoscópio descartáveis	53
GR	Αναλώσιμες ενδοσκοπικές βελόνες έγχυσης	55
TR	Tek Kullanımlık Endoskop Enjeksiyon İğneleri	57

EXPLANATION OF SYMBOLS USED ON LABELS AND INSTRUCTIONS FOR USE

Symbol	Explanation
	Single use only
	Date of manufacture
	Manufacturer
	Use-by-date
	Catalogue number
	Batch code
	CE mark
	Keep dry
	Caution
	Consult instructions for use
	Sterilized by Ethylene Oxide
	Do Not Re-Sterilize
	Authorized representative in the European Community
	Authorized representative in the United Kingdom
	Do Not Use if packaging or product damage is evident
	Diameter
	Applied to channel diameter
	Working length
	Pack Quantity
	Medical Device

INTENDED USE

The device enters into the human body alimentary canal through the endoscope to perform endoscopic vascular or submucosal injection in the GI tract.

CAUTION

Flexible endoscopy procedures should only be performed by persons having adequate training and familiarity with endoscopic techniques. Consult the medical literature relative to techniques, contraindications, complications and hazards prior to any endoscopic procedure.

- The ECC^{INJECT} are single use, disposable medical devices; they are not intended for reuse. Any institution, practitioner or third party who reprocesses, refurbishes, remanufactures, sterilizes, and/or reuses these disposable devices must bear full responsibility for their safety and effectiveness.
- The manufacturer did not design these devices to be reprocessed or reused, and therefore cannot verify that reprocessing can clean and/or sterilize or maintain the structural integrity of the device to ensure patient and/or user safety.
- Do not attempt to actuate the ECC^{INJECT} with the catheter in a straightened position. The ECC^{INJECT} has been engineered and manufactured with a specific "preload" that is necessary to assure the full projection and full retraction of the needle from the distal hub when the device is in very tortuous configurations. If the ECC^{INJECT} is actuated with the catheter in a straightened position, the "preload" may cause the needle's inner catheter to collapse and not allow for adequate injection. When handling the ECC^{INJECT}, always be aware of the distal tip of the needle, and visually confirm that the needle is fully retracted into the device.

Do not attempt to activate or test the ECC^{INJECT} outside of the endoscope if the device is coiled more than once. It is possible that this may cause the ECC^{INJECT} to become lodged proximal to the distal tip spring.

- Compatibility with medicinal substances:
 - All of the ECC^{INJECT} included in this IFU (Instructions for Use) are compatible with the following medicinal substances and solutions: Sodium Morrhuate, Sodium Tetradecyl Sulfate, Ethynolamine Oleate, Epinephrine, Botulinum Toxin Type A, Polidocanol, Ethanol, Hypertonic Saline and Staining Media.
 - Note: Sodium Morrhuate, Sodium Tetradecyl Sulfate, Polidocanol, and Ethanol have not been approved by the United States Federal Food and Drug Administration and should not be used in endoscopic procedures.
 - The compatibility of the ECC^{INJECT} with other media has not been verified and therefore should not be attempted.

CONTRAINDICATIONS

Contraindications include those specific to any endoscopic procedure, as well as those specific to endoscopic needle therapy.

DIRECTION OF USE

The operator of the instrument must be a physician or medical personnel under the supervision of a physician and must have received sufficient training in clinical endoscopic technique. This manual, therefore, does not explain or discuss clinical endoscopic procedures. It only describes basic operation and precautions related to the operation of this instrument.

WARNING

- When using the instrument, always wear appropriate personal protective equipment. Otherwise, blood, mucous and other potentially infectious material from the patient could pose an infection control risk, appropriate personal protective equipment may include; Eye wear, a face mask, moisture-resistant clothing and chemical-resistant gloves that fit properly and are long enough so that your skin is not exposed.
- Do not insert the instrument into the endoscope unless you have a clear endoscopic field of view, if you cannot see the distal end of the insertion portion in the endoscopic field of view, do not use it. This could cause patient injury, such

CONTACT:

Manufacturer:
Yangzhou Fartley Medical Instrument
Technology Co., Ltd. - Beizhou Road,
Lidian Town, Guangling District,
Yangzhou, Jiangsu 225106, China

EC REP
Prolinx GmbH - Brehmstr.56, D-40239
Duesseldorf, Germany
SRN: DE-AR-000005129

UK REP
Prolinx Global LTD
27 Old Gloucester Street, London
WC1N 4AX - UK - MHRA
Account Number: 0000013089

International distributor:

mdd Company GmbH - Switzerland
Rengglochstr. 38 CH-6012 Obernau
E-Mail: info@mddcompany-gmbh.com
www.mddcompany-gmbh.com

as perforation, bleeding, or mucous membrane damage, it may also damage the endoscope and/or instrument.

- Do not angulate the endoscope's bending section or operate the forceps elevator abruptly while the distal end of the insertion portion is extended from the distal end of the endoscope, this could cause patient injury, such as perforation, bleeding, or mucous membrane damage.
- Do not force the distal end of the insertion portion against body cavity tissue. This could cause patient injury, such as perforation, bleeding, or mucous membrane damage.

Insertion into the Endoscope

WARNING

- Do not insert the instrument into the endoscope when the needle is extended.
- The distal end of the insertion portion may extend from the endoscope tip abruptly, this could cause patient injury, such as perforation, bleeding, or mucous membrane damage. It could also damage the endoscope and/or instrument.
- When inserting the instrument into the endoscope, hold the holder portion and slider firmly. Otherwise, the needle may extend from the distal end of endoscope abruptly. This could cause patient injury, such as perforation, bleeding, or mucous membrane damage, it could also damage the endoscope or instrument.
- Do not force the instrument if resistance to insertion is encountered. Reduce the angulation or lower the forceps elevator of the endoscope until the instrument passes smoothly. This could cause patient injury, such as perforation, bleeding, or mucous membrane damage, it could also damage the endoscope and/or instrument.
- Do not advance or extend the instrument abruptly, this could cause patient injury, such as perforation, bleeding, or mucous membrane damage, it could also damage the endoscope and/or instrument.


COMPATIBILITY

Our mission is to produce safe therapeutic endoscopic accessories with the highest level of quality possible.

All Medical accessories have been tested in conjunction with flexible medical endoscopes for function and compatibility. The voice to the customer is important to keep consistent quality. We are in close contact with customers which confirms the stable quality and compatibility. The minimum requirement for using our ECC[®] with flexible Endoscope is to guarantee that biopsy channel / operation channel is minimum 0,2 mm bigger than the outer diameter of the ECC[®] sheath.

CAUTION

- When inserting the instrument into the endoscope. Hold it close to the biopsy valve and keep it as straight as possible relative to the biopsy valve. Otherwise the instrument could be damaged.
 - Insert the instrument slowly, abrupt insertion could damage the endoscope and/or instrument.
1. If the endoscope is equipped with a forceps elevator, lower the forceps elevator.

2. Pull the slider until it stops to retract the needle into the distal end of the tube.

3. With the needle retracted, carefully insert the instrument into the biopsy valve of the endoscope.

4. Advance the instrument until the distal end of the insertion portion appears within the endoscopic field of view.

Piercing tissue and injecting the medicine

WARNING

- Do not push the slider abruptly, otherwise the needle will be rapidly extended from the distal end of the tube, this could result in patient injury, such as perforation, bleeding, or mucous membrane damage. It could also damage the instrument.
- While piercing tissue and injecting medication, always hold the holder portion and slider, otherwise, this could result in patient injury, such as hemorrhages.
- Do not pierce tissue with excessive force. This could result in patient injury, such as perforation, bleeding, or mucous membrane damage.
- When the slider clicks into position after being pushed, keep the insertion portion and the bending section of the endoscope as straight as possible. Otherwise, the tube may buckle which could cause in to pierce tissue and inject medication.
- Push the slider into the holder portion until it clicks into position to extend the needle from the distal end of the tube.
- Connect a syringe filled with a medicine to the injection port of the slider.
- Slightly depress the syringe's plunger and confirm that medicine comes out of the distal end of the needle.
- Pierce the target tissue with the needle. Depress the plunger to inject the medicine.
- Withdraw the needle from the target tissue.
- Pull the slider until it stops to retract the needle into the distal end of the tube.


Withdrawing the instrument from the endoscope

WARNING: Do not withdraw the instrument from the endoscope quickly. This could scatter blood, mucous or other patient debris and pose an infection control risk.

CAUTION: Do not withdraw the instrument from the endoscope while the needle is extended. This could damage the endoscope or instrument.

If the endoscope is equipped with a forceps elevator, lower the forceps elevator.

Withdraw the instrument from the endoscope.

CE
0197

Warning: An issued or revision date for these instructions is included for the user's information. In the event that two years have elapsed between this date and product use, the user should contact Company Fartley to determine if additional information is available. Serious incidents that have occurred in relation to this medical device should be reported to the manufacturer and competent authority in the country where incident occurred.

Product Disposal: After use, this product may be a potential biohazard. Handle and dispose of in accordance with accepted medical practice and applicable local, state and federal laws and regulations.

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The scope of delivery of the medical device Disposable Endoscope Injection Needle includes one of the MD variations at the customer's choice in accordance with the Name:


1. Disposable Endoscope Injection Needle **IN-O18-19-4-180**, containing:
 - 1.1 Disposable Endoscope Injection Needle **IN-O18-19-4-180**, 10 pcs/pack
 - 1.2 Instruction manual, 1 pc.
2. Disposable Endoscope Injection Needle **IN-O18-19-5-180**, containing:
 - 2.1 Disposable Endoscope Injection Needle **IN-O18-19-5-180**, 10 pcs/pack
 - 2.2 Instruction manual, 1 pc.
3. Disposable Endoscope Injection Needle **IN-O18-19-6-180**, containing:
 - 3.1 Disposable Endoscope Injection Needle **IN-O18-19-6-180**, 10 pcs/pack
 - 3.2 Instruction manual, 1 pc.
4. Disposable Endoscope Injection Needle **IN-O18-22-4-180**, containing:
 - 4.1 Disposable Endoscope Injection Needle **IN-O18-22-4-180**, 10 pcs/pack
 - 4.2 Instruction manual, 1 pc.
5. Disposable Endoscope Injection Needle **IN-O18-22-5-180**, containing:
 - 5.1 Disposable Endoscope Injection Needle **IN-O18-22-5-180**, 10 pcs/pack
 - 5.2 Instruction manual, 1 pc.
6. Disposable Endoscope Injection Needle **IN-O18-22-6-180**, containing:
 - 6.1 Disposable Endoscope Injection Needle **IN-O18-22-6-180**, 10 pcs/pack
 - 6.2 Instruction manual, 1 pc.
7. Disposable Endoscope Injection Needle **IN-O18-23-4-180**, containing:
 - 7.1 Disposable Endoscope Injection Needle **IN-O18-23-4-180**, 10 pcs/pack
 - 7.2 Instruction manual, 1 pc.
8. Disposable Endoscope Injection Needle **IN-O18-23-5-180**, containing:
 - 8.1 Disposable Endoscope Injection Needle **IN-O18-23-5-180**, 10 pcs/pack
 - 8.2 Instruction manual, 1 pc.
9. Disposable Endoscope Injection Needle **IN-O18-23-6-180**, containing:

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
























- Individual packaging label may be one of two available versions, depending on the date of manufacture of the medical device.

The label version of the MD was manufactured before 01/01/2023 contains the following symbols/information about the medical device:

- sterilization control indicator;
- inscription in English “E.O.GAS BROWN-PROCESSED IN TO GAS”, meaning: “brown indicator color – processed by ethylene oxide gas”;
- inscription in English: «Do not use if pack is damaged»;
- brand name of the medical device: **ECO_{med}**;
- information about the company – international distributor: name, address, website;
- catalogue number of the variation, symbol «catalogue number»;
- number of the medical device according to the internal manufacturer’s catalogue;
- batch code, symbol «batch code»;
- date of manufacture in format: yyyy-mm-dd;
- symbol «date of manufacture»;
- use-by-date in format: yyyy-mm-dd;
- symbol «use-by-date»;
- number of pcs medical device in the package, symbol «pack quantity»;
- schematic picture of the MD with color coding depending on the needle diameter (purple for 19G, green for 22G, blue for 23G, pink for 25G);
- needle diameter in mm in format: Ø X.X mm;
- needle length in mm in format: Length: X mm;
- working length in cm in format: XXX cm;

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The label version of the MD was manufactured before 01/01/2023 contains the following symbols/information about the medical device:

- sterilization control indicator;
- inscription in English “E.O.GAS BROWN-PROCESSED IN TO GAS”, meaning: “brown indicator color – processed by ethylene oxide gas”;
- inscription in English: «Do not use if pack is damaged»;
- brand name of the medical device: **ECCO^{med}**;
- information about the company – international distributor: name, address, website;
- catalogue number of the variation, symbol «catalogue number»;
- number of the medical device according to the internal manufacturer’s catalogue;
- batch code, symbol «batch code»;
- date of manufacture in format: yyyy-mm-dd;
- symbol «date of manufacture»;
- use-by-date in format: yyyy-mm-dd;
- symbol «use-by-date»;
- number of pcs medical device in the package, symbol «pack quantity»;
- schematic picture of the MD with color coding depending on the needle diameter (purple for 19G, green for 22G, blue for 23G, pink for 25G);
- needle diameter in mm in format: Ø X.X mm;
- needle length in mm in format: Length: X mm;
- working length in cm in format: XXX cm;

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- Consumer packaging label may be one of two available versions, depending on the date of manufacture of the medical device.

The label version of the MD was manufactured before 01/01/2023 contains the following symbols/information about the medical device:

- brand name of the medical device: ;
- needle diameter in gage in the format: XX Gauge;
- information about the company – international distributor: name, address, website;
- catalogue number, symbol «catalogue number»;
- number of the medical device according to the manufacturer's catalogue;
- batch code, symbol «batch code»;
- date of manufacture in format: yyyy-mm-dd;
- symbol «date of manufacture»;
- use-by-date in format: yyyy-mm-dd;
- symbol «use-by-date»;
- number of medical devices per box, symbol «pack quantity»;
- schematic picture of the MD with color coding depending on the needle diameter (purple for 19G, green for 22G, blue for 23G, pink for 25G);
- needle diameter in mm in format: Ø X.X mm;
- needle length in mm in format: Length: X mm;

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