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To whom it may concern

STATEMENT

To be submitted to the Competent Health Authorities in Russian Federation for the medical device «HILZO® gastroenterological stents in various versions».

We, BCM Co., Ltd, 178 Deoksan-ro, Ilsanseo-gu, Goyang-si, Gyeonggi-do, Republic of Korea hereby state that attached document «Operational documentation »

are true and identical copies of the original documents.

DATE: 21 / 11 / 2023

NAME: Myung, Byung-cheol

POSITION: CEO

Stamp and signature

BCM Co., Ltd.
178, Deoksan-ro, Ilsanseo-gu,
Goyang-si, Gyeonggi-do, KOREA
Tel +82-31-926-2727
CEO Myung, Byung-Cheol

We hereby certify that the signature appearing herein has been duly registered at this chamber

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
LEE, JAEHOON -----

attorney-in-fact of
BCM CO., LTD.
CEO MYUNG, BYUNG-CHEOL -----

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appeared before me and admitted
said principal`s subscription to
the attached STATEMENT. -----

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This is hereby attested on this
24th day of Nov. 2023 at this office.

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(Cheongjin-dong, Samsong Bldg.) 701, 702,
34, Jong-ro 3-gil, Jongno-gu, Seoul,
Republic of Korea

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Signature of the Notary Public
Lee Kyung Taek

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This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public Since
7, Feb. 2020 Under Law No.243.

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This public document

2. has been signed by Lee Kyung Taek

3. acting in the capacity of Notary Public

4. bears the seal/stamp of Law Firm Idea AND NOTARY OFFICE INC.

Certified

To verify the Apostille, please refer to the website below.
<https://www.apostille.go.kr>

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
ENGLISH

SYMBOLS

CE mark is a mandatory conformity marking for certain products sold within the European Economic Area (EEA).

Authorized representative in the European Community

Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.

Indicates a medical device that should not be used if the package has been damaged or opened.

Indicates the date after which the medical device is not to be used.

Indicates the date when the medical device was manufactured.

Indicates the manufacturer's catalogue number so that the medical device can be identified.

Indicates a medical device that has been sterilized using ethylene oxide.

Indicates a medical device that is not to be re-sterilized.

Indicates the need for the user to consult the instructions for use for important cautionary information.

Indicates the temperature limits to which the medical device can be safely exposed

Indicates the medical device manufacturer, as defined in EU Directives 90/385/EEC, 93/42/EEC and 98/79/EC.

Indicates the manufacturer's batch code so that the batch or lot can be identified.

Indicates the manufacturer's serial number so that a specific medical device can be identified.

MR conditional

1. Description


The HILZO® Biliary Stent system consists of a delivery system preloaded with a self-expanding biliary metal stent. The stent is made of Nitinol wire knitted in a tubular mesh configuration. This design configuration makes a stent flexible and self-expandable. The stent has 6 radiopaque markers.

Table 1. Stent Model

Model Name
HILZO® Biliary Uncovered Stent
HILZO® Biliary Covered Stent

HILZO® Full Covered Biliary Stents can be removed; (see Warnings).

Figure 2. HILZO® Biliary Stent System

The stent is mounted on an inner catheter and is constrained by an outer tube. The outer tube is pulled back by immobilizing the metal tube holder with one hand, grasping the handle with the other hand, and gently sliding the handle along the metal tube towards the metal tube holder. With a retraction of the outer tube, the stent is released. When the Stent is deployed, the stent self-expands and imparts an outward radial force on the luminal surface of biliary to establish patency. For biliary stenting, two different approaches are utilized: Percutaneous approach and Endoscopic approach.

HILZO® Stent system has two different delivery systems depending on the approach.

Table 2. Recommended delivery system by approach

Approach	Usable Length	Delivery system type
Percutaneous approach	50cm	Percutaneous
Endoscopic approach	180cm	Endoscopic

2. Indication for Use

The HILZO® Biliary Stent is intended for use in the palliative treatment of biliary strictures produced by malignant neoplasm. Fully covered Biliary Stent is intended for the treatment of benign biliary strictures.

WARRANTY

BCM Co., LTD. warrants that reasonable care has been applied within the design and subsequent manufacturing process of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether expressed or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond BCM's control directly affect the instrument and the results obtained from its use. BCM's obligation under this warranty is limited to the replacement of this instrument and BCM shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. BCM neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. BCM assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, expressed or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instruments.

3. Contraindication

The HILZO® Biliary Stent is contraindicated for, but are not limited to:

- Patients with coagulopathy
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bared Stents is contraindicated (see Warnings).
- Suspected or impending perforation.
- Patients for whom endoscopic techniques are contraindicated.
- Biliary obstruction preventing either endoscopic or percutaneous cholangiography.

4. Warnings

- The safety and efficacy of this device for use in the vascular system has not been established.
- The device should be used with caution and only after careful consideration in patients

with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.

- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration or/and fracture.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the delivery system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- HILZO® Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- HILZO® Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Uncovered Stents should not be removed once fully deployed; see Contraindications.
- Visually inspect the system for any sign of damage. **DO NOT USE** if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

5. Potential complications

The following complications have been reported in the literature for biliary prosthesis. Potential complications associated with the use of HILZO® Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Intestinal perforation

Post Stent Placement Complications

- Bleedings
- Pain
- Nausea
- Vomiting
- Inflammation
- Mucosal hyperplasia
- Recurrent obstructive jaundice
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Tumor in-growth through stent
- Tumor over-growth around ends of stent
- Fever
- Ulceration
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Pancreatitis
- Cholangitis/Cholestasis/Cholecystitis
- Constipation
- Diarrhea
- Liver abscess
- Peritonitis
- Sludge occlusion

6. Required equipment

Percutaneous Placement

- 0.035" (0.89mm) guidewire at least 180cm long (preferably stiff or extra stiff)
- Introducer sheath appropriately sized for HILZO® stent and delivery system

Endoscopic Placement

- 0.035" (0.89mm) guidewire at least 450cm long
- Endoscope with minimum 3.2mm working channel

7. Precautions

Read the entire User's Manual thoroughly before using this device. It should only be used by or under the supervision of physicians thoroughly trained in the placement of stents. A thorough understanding of the techniques, principles, clinical applications and risks associated with this procedure is necessary before using the device.

- Care should be taken when removing the delivery system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.

- The packaging and the device should be inspected prior to use.
- Do not attempt to reload deployed stents into the delivery system.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The HILZO® Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The HILZO® Stent is intended for single use only. Do not resterilize and/or reuse the device.
- If the stent is not fully deployed or delivery system is not operated as intended by this manual, stop the operating and contact BCM.
- Be careful not to contaminate the product during package open.

B. Pre-Procedure Note

- Carefully remove the delivery system from the packaging.
- Visually inspect the device for damage or defects.

3. Procedure

- **Delivery system with a 50cm usable length for percutaneous approach.**
Percutaneous transhepatic cholangiography (PTC) should be performed prior to the placement of the HILZO® stent to examine the biliary tract morphology and extent of the stricture.
 - **Delivery system with a 180cm usable length for Endoscopic approach.**
Endoscopic retrograde cholangiopancreatography (ERCP) should be performed prior to the placement of the HILZO® Stent to examine the biliary tract morphology and extent of the stricture.
- ① Examine stricture endoscopically and fluoroscopically
 - a) Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
 - b) The internal luminal diameter should be measured exactly with an endoscope and/or a fluoroscope.
 - ② Stent Size selection
 - a) Measure the length of the target stricture.
 - b) Select a stent size that is 20mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
 - c) Check stent size and specification on label.
 - ③ Stent Deployment Preparation
 - The HILZO® Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass a 0.035" (0.89 mm) guidewire to the level of the stricture.
 - A. Fluoroscopy Procedure
 - a) Under fluoroscopic/endoscopic guidance, insert a guide wire across the stricture to where the stent delivery system will be placed over the guide wire.
 - b) Ensure that the outer tube lock system is locked to prevent premature stent deployment. The outer tube lock system is locked by rotation in a clockwise direction
 - B. Endoscopy Procedure
 - a) Under the fluoroscopic/endoscopic guidance, insert an endoscope to the level of the obstruction, then introduce the guide wire through the working channel of the endoscope. Advance until the guide wire passes the target stricture to where the stent delivery system will be placed over the guide wire.
 - b) Ensure that the outer tube lock system is locked to prevent premature stent deployment. The outer tube lock system is locked by rotation in a clockwise direction
 - ④ Stent Deployment Procedure


Figure 3

PRECAUTION: Do not twist delivery system or employ a boring motion during the deployment as this may affect positioning and ultimate function of stent.

- a) Under fluoroscopic and/or endoscopic guidance, position the delivery system to the center of the target stricture exactly.
- b) Once the delivery system is in the correct position for deployment, unlock the outer tube lock system by turning it more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the metal tube holder with one hand and grasp the handle with the other hand. Gently slide the handle back along the metal tube towards the metal tube holder.

- d) When the center X-ray marker reaches the center of target stricture, continue pulling back the handle until the stent is fully deployed. (See figure 3, 4)


Figure 4

CAUTION: Do not push forward or pull backward the metal tube holder with the stent partially deployed. The metal tube holder must be securely immobilized. Inadvertent movement of the metal tube holder may cause misalignment of the stent and possible damage to the bile duct.

⑤ After Stent Deployment

- Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- Carefully remove the delivery system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3-5 minutes to allow further stent expansion.

10. Perform routine post implant procedures

- Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- Observe the patient for development of any complication.

11. Instructions for removal of HILZO® Covered Stents (see Warnings)

Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

To reposition a HILZO® Covered Stent immediately after deployment, use forceps or a snare to grasp the end of the stent and gently adjust to the correct placement.

Please note: the stent can only be repositioned and/or removed proximally.

Reuse Precaution Statement

Contents supplied STERILE (ethylene oxide (EO)). Do not use if sterile barrier is damaged. In the event of damaged packaging, call your BCM Co., Ltd. representative. For single patient use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross infection, including, but not limited to, the transmission of infectious diseases from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

Storage: Store at room temperature.

Disposal Requirements: After use, dispose product and packaging in accordance with hospital, administrative and/or local government policy for biohazardous medical waste.

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MRI Conditional

1. Description

The HILZO® Esophageal Stent system consists of a delivery system preloaded with a self-expanding esophageal metal stent. The stent is made of Nitinol wire knitted in a tubular mesh configuration. This design configuration makes a stent flexible and self-expandable. The stent has 6 radiopaque markers.

The stent is mounted on an inner catheter and is constrained by an outer tube. The outer tube is pulled back by immobilizing the metal tube holder with one hand, grasping the handle with the other hand, and gently sliding the handle along the metal tube (2nd inner catheter) towards the metal tube holder. With a retraction of the outer sheath, the stent is released. When the Stent is deployed, the stent self-expands and imparts an outward radial force on the luminal surface of the esophageal to establish patency.

Figure 2. HILZO® Esophageal Delivery System

For esophageal stenting, two different approaches are utilized: OTW (Over-The-Wire) and TTS (Through-The-endoscope) approach.

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Table 2. Delivery system Type

Delivery system type	Usable Length	Approach
OTW (Over-the-Wire)	70cm	The stent delivery system is passed over the guidewire into the esophagus
TTS (Through The endoscope)	180cm	The stent delivery system is inserted into the esophagus through an endoscope

2. Indication for Use

The HILZO® Esophageal Stent is intended for maintaining esophageal luminal patency in esophageal stricture caused by intrinsic and/or extrinsic malignant tumors, occlusion of concurrent esophageal fistulas, refractory benign esophageal strictures.

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3. Contraindication

The HILZO® Esophageal Stent is contraindicated for, but is not limited to:

- Placement in polypoid lesions
- Patient with bleeding disorder
- Strictures that do not allow passage of a guidewire
- Any use other than those specifically outlined under indications for use.
- Suspected or impending perforation
- Strictures that cannot be dilated enough to pass the endoscope or the delivery system
- Placement of the proximal end of the stent within 2 cm of the cricopharyngeal muscle

4. Warnings

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, significant preexisting pulmonary or cardiac disease.

- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration and/or fracture.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the delivery system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- HILZO® Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- HILZO® Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Uncovered/bare Stents should not be removed once fully deployed; see Contraindications.
- Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

5. Potential complications

Potential complications associated with the use of HILZO® Stent may include, but are not limited to:

- Bleeding
- Pain
- Perforation
- Stent misplacement or migration
- Stent occlusion
- Stent occlusion due to tumor in-growth through Stent
- Stent occlusion due to tumor over-growth around ends of Stent
- Tumor in-growth
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Intestinal obstruction secondary to migration/Haematemesis
- Airway Compressions
- Reflux
- Food bolus impaction
- Esophagitis
- Recurrent Dysphagia
- Fistula formation
- Ulcerations
- Aspiration
- Edema
- Nausea

6. Required equipment

- Fluoroscope and/or endoscope with minimum 3.7mm working channel
- 0.035" (0.89 mm) Guidewire

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- Care should be taken when removing the delivery system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Do not attempt to reload deployed stents into the delivery system.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The HILZO® Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The HILZO® Stent is intended for single use only. Do not resterilize and/or reuse the device.
- Be careful not to contaminate the product during package open.

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9. Procedure

- ① Examine stricture endoscopically and fluoroscopically
 - a) Carefully examine both the proximal and distal segment of the stricture endoscopically and/or fluoroscopically.
 - b) The internal luminal diameter should be measured exactly with an endoscope and/or a fluoroscope.
- ② Stent Size selection
 - a) Measure the length of the target stricture.
 - b) Select a stent size that is 20mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
 - c) Check stent size and specification on label.
- ③ Stent Deployment Preparation
 - The HILZO® Esophageal Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass guidewire to the level of the stricture.
 - a) Under fluoroscopic/endoscopic guidance, insert a guide wire across the stricture to where the stent delivery system will be placed over the guide wire.
 - b) Ensure that the outer tube lock system is locked to prevent premature stent deployment. The outer tube lock system is locked by rotation in a clockwise direction.
- ④ Stent Deployment Procedure


Figure 3

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- a) Under fluoroscopic and/or endoscopic guidance, position the delivery system to the center of the target stricture exactly.
- b) Once the delivery system is in the correct position for deployment, unlock the outer tube lock system by turning it more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the metal tube holder with one hand and grasp the handle with the other hand. Gently slide the handle back along the metal tube towards the metal tube holder.
- d) When the center X-ray marker reaches the center of target stricture, continue pulling back the handle until the stent is fully deployed. (See Figure 3, 4)


Figure 4

CAUTION: Do not push forward or pull backward the metal tube holder with the stent partially deployed. The metal tube holder must be securely immobilized. Inadvertent movement of the metal tube holder may cause misalignment of the stent and possible damage to the esophagus.

- ⑤ After Stent Deployment
 - a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
 - b) Carefully remove the delivery system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion

c) Balloon dilatation inside the stent can be performed on demand.

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- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After Implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complication.

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Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent. To reposition a HILZO® Covered Stent immediately after deployment, use forceps or a snare to grasp the end of the stent and gently adjust to the correct placement. Please note: the stent can only be repositioned and/or removed proximally.

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Storage: Store at room temperature.

Disposal Requirements: Upon completion of procedure, device must be properly disposed in accordance with local or hospital regulations for biohazardous medical waste.

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MRI Conditional

1. Description

The HILZO® Pyloric/Duodenal Stent system consists of a delivery system preloaded with a self-expanding Pyloric/Duodenal metal stent. The stent is made of Nitinol wire knitted in a tubular mesh configuration. This design configuration makes a stent flexible and self-expandable. The stent has 6 radiopaque markers.

Table 1. Stent Model

Model Name
HILZO® Pyloric/Duodenal Uncovered Stent
HILZO® Pyloric/Duodenal Covered Stent

Figure 2. Delivery system

The stent is mounted on an inner catheter and is constrained by an outer catheter. The outer sheath is pulled back by immobilizing the metal tube holder with one hand, grasping the handle with the other hand, and gently sliding the handle along the metal tube towards the metal tube holder. With a retraction of the outer tube, the stent is released. When the Stent is deployed, the stent self-expands and imparts an outward radial force on the luminal surface of pyloric/duodenum to establish patency. The delivery system may be utilised for endoscopic delivery as the profile of the delivery system is compatible with a minimum 3.7mm working channel.

2. Indication for Use

The HILZO® Pyloric/Duodenal is intended for the palliative treatment of gastroduodenal obstruction produced by malignant strictures.

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3. Contraindication

The HILZO® Pyloric/Duodenal Stent is contraindicated for, but is not limited to:

- Intra-abdominal abscess.
- Patients with coagulopathy
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under indications for use.
- Removal or repositioning of fully deployed uncovered/bared Stents is contraindicated (see Warnings).
- Suspected or impending perforation.
- Enteral ischemia

4. Warnings

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis.
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration or/and fracture.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the delivery system to organic solvent (e.g. Alcohol)

- Do not use with Ethiodol or Lipiodol contrast media.
- HILZO® Covered Stents can be removed where the treating doctor's clinical assessment of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.
- HILZO® Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Uncovered/bare Stents should not be removed once fully deployed; see Contraindications.
- Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

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Past Stent Placement Complications

- Bleeding
- Pain
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Stent occlusion due to tumor in-growth through stent
- Stent occlusion due to tumor over-growth around ends of stent
- Tumor in-growth
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Infection
- Peritonitis
- Sludge occlusion
- Ulceration
- Jejunal stenosis

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- Care should be taken when removing the delivery system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or stent migration.
- The packaging and the device should be inspected prior to use.
- Do not attempt to reload deployed stents into the delivery system.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The HILZO® Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The HILZO® Stent is intended for single use only. Do not resterilize and/or reuse the device.

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- b) Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
 - b) The internal luminal diameter should be measured exactly with an endoscope and/or a fluoroscope.
- ② Stent Size selection
- a) Measure the length of the target stricture.
 - b) Select a stent size that is 20 mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
- ③ Stent Deployment Preparation
- The HILZO® Pyloric/Duodenal Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass a guidewire to the level of the stricture.
 - a) Under fluoroscopic guidance, insert a guide wire across the stricture to where the stent delivery system will be placed over the guide wire.
 - b) Ensure that the outer tube lock system is locked to prevent premature stent deployment. The outer tube lock system is locked by rotation in a clockwise direction.
 - c) Flush the inner lumen of delivery system.
- ④ Stent Deployment Procedure


Figure 3

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- a) Under fluoroscopic and/or endoscopic guidance, position the delivery system to the center of the target stricture exactly.
- b) Once the delivery system is in the correct position for deployment, unlock the outer tube lock system by turning it more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the metal tube holder with one hand and grasp the handle with the other hand. Gently slide the handle back along the metal tube towards the metal tube holder.
- d) When the center X-ray marker reaches the center of target stricture, continue pulling back the handle until the stent is fully deployed. (See figure 3, 4)


Figure 4

CAUTION: Do not push forward or pull backward the metal tube holder with the stent partially deployed. The metal tube holder must be securely immobilized. Inadvertent movement of the metal tube holder may cause misalignment of the stent and possible damage to pylorus or duodenum.

⑨ After Stent Deployment

- a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
- b) Carefully remove the delivery system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion
- c) Balloon dilatation inside the stent can be performed on demand.

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- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.
- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complication.

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To reposition a HILZO® Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement. Please note: the stent can only be repositioned and/or removed proximally.

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MR conditional


1. Description

The HILZO® Colonic Stent system consists of a delivery system preloaded with a self-expanding metal stent. The stent is made of Nitinol wire knitted in a tubular mesh configuration. This design configuration makes a stent flexible and self-expandable. The stent has 6 radiopaque markers.

Table 1. Stent Model

Model Name
HILZO® Colonic Uncovered Stent
HILZO® Colonic Covered Stent

Figure 2. Delivery system

The stent is mounted on an inner catheter and is constrained by an outer tube. The outer tube is pulled back by immobilizing the metal tube holder with one hand, grasping the handle with the other hand, and gently sliding the handle along the metal tube towards the metal tube holder. With a retraction of the outer tube, the stent is released. When the Stent is deployed, the stent self-expands and imparts an outward radial force on the luminal surface of the colon to establish patency. The delivery system may be utilised for endoscopic approach as the profile of the delivery system is compatible with a minimum 3.7mm working channel.

2. Indication for Use

The HILZO® Colonic Stent is intended for the palliative treatment of colonic strictures caused by malignant neoplasm and to relieve large bowel obstruction.

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3. Contraindication

The HILZO® Colonic Stent is contraindicated for, but is not limited to:

- Intra-abdominal abscess.
- Patients with coagulopathy
- Strictures that do not allow passage of a guidewire.
- Any use other than those specifically outlined under Indications for use.
- Removal or repositioning of fully deployed uncovered/bared Stents is contraindicated (see Warnings).
- Suspected or impending perforation.
- Enteral ischemia

4. Warnings

- The device should be used with caution and only after careful consideration in patients with elevated bleeding times, coagulopathies, or in patients with radiation colitis or proctitis
- Chemoradiation therapy or radiotherapy alone may lead to tumor shrinkage and subsequent stent migration or/and fracture.
- The stent contains nickel, which may cause an allergic reaction in individuals with nickel sensitivity.
- Do not expose the delivery system to organic solvent (e.g. Alcohol)
- Do not use with Ethiodol or Lipiodol contrast media.
- HILZO® Covered Stents can be removed where the treating doctor's clinical assessment

of the stented stricture indicates the Stent can be safely removed. Caution should be exercised in deciding to and when removing the Stent.

- HILZO® Covered Stent cannot be removed when there is tumor in-growth/over-growth/occlusion of the Stent lumen.
- Uncovered/bare Stents should not be removed once fully deployed; see Contraindications.
- Visually inspect the system for any sign of damage. DO NOT USE if the system has any visible signs of damage. Failure to observe this precaution may result in patient injury.

5. Potential complications

Potential complications associated with the use of HILZO Stent may include, but are not limited to:

Procedural Complications

- Bleeding
- Stent misplacement or inadequate expansion
- Pain
- Death (other than that due to normal disease progression)
- Intestinal Perforation

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- Bleeding
- Pain
- Perforation
- Bowel impaction
- Stent misplacement or migration
- Stent occlusion
- Stent occlusion due to tumor in-growth through stent
- Stent occlusion due to tumor over-growth around ends of stent
- Tumor ingrowth
- Fever
- Foreign body sensation
- Death (other than that due to normal disease progression)
- Sepsis
- Constipation
- Diarrhea
- Infection
- Peritonitis
- Sludge occlusion
- Ulceration
- Symptoms of tenesmus or urgency/incontinence

6. Required equipment

- Fluoroscope and/or endoscope (with a channel size of 3.7 mm or larger)
- 0.035" (0.89 mm) Guidewire

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- Care should be taken when removing the delivery system and guidewire immediately after stent deployment since this may result in stent dislodgement if the stent has not been adequately deployed.
- Care should be taken when performing dilation after the Stent has been deployed as this may result in perforation, bleeding, Stent dislodgement or Stent migration.
- The packaging and the device should be inspected prior to use.
- Do not attempt to reload deployed stents into the delivery system.
- Use of fluoroscopy is recommended to ensure correct placement of the device.
- Check the expiration date "Use by". Do not use the device beyond the use by date.
- The HILZO® Stent is supplied sterile. Do not use if the packaging is opened or damaged.
- The HILZO® Stent is intended for single use only. Do not resterilize and/or reuse the device.
- Be careful not to contaminate the product during package open

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9. Procedure

- ① Examine stricture endoscopically and fluoroscopically

- a) Carefully examine both the proximal and distal segment of stricture endoscopically and/or fluoroscopically.
 - b) The internal luminal diameter should be measured exactly with an endoscope and/or a fluoroscope.
- ② Stent Size Selection
- a) Measure the length of the target stricture.
 - b) Select a stent size that is 20 mm longer than the measured length of the stricture in order to cover fully both ends of the lesion.
 - c) Check stent size and specification on label.
- ③ Stent Deployment Preparation
- The HILZO® colonic Stent can be placed with the aid of fluoroscopy and/or endoscopy.
 - Pass a guidewire to the level of the stricture.
- a) Under fluoroscopic/Endoscope guidance, insert a guide wire across the stricture to where the stent delivery system will be placed over the guide wire.
 - b) Ensure that the outer tube lock system is locked to prevent premature stent deployment. The outer tube lock system is locked by rotation in a clockwise direction.
- ④ Stent Deployment Procedure


Figure 3

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- a) Under fluoroscopic and/or endoscopic guidance, position the delivery system to the center of the target stricture exactly.
- b) Once the delivery system is in the correct position for deployment, unlock the outer tube lock system by turning it more than twice in an anti-clockwise direction.
- c) To begin stent deployment, immobilize the metal tube holder with one hand and grasp the handle with the other hand. Gently slide the handle back along the metal tube towards the metal tube holder.
- d) When the center X-ray marker reaches the center of the target stricture, continue pulling back the handle until the stent is fully deployed. (See figure 3, 4)


Figure 4

CAUTION Do not push forward or pull backward the metal tube holder with the stent partially deployed. The metal tube holder must be securely immobilized. Inadvertent movement of the metal tube holder may cause misalignment of the stent and possible damage to the colon.

- ⑤ After Stent Deployment
- a) Examine the stent fluoroscopically and/or endoscopically to confirm expansion.
 - b) Carefully remove the delivery system, guidewire and endoscope from the patient. If excessive resistance is felt during removal, wait 3~5 minutes to allow further stent expansion.
 - c) Balloon dilatation inside the stent can be performed on demand.

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- a) Assess the size and stricture of the Stent lumen. A Stent may require up to 1 to 3 days to expand fully.

- b) Doctor's experience and discretion can determine the appropriate drug regimen for each patient.
- c) After implantation, patient should remain on a soft diet until otherwise determined by the treating doctor.
- d) Observe the patient for development of any complication.

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Visually examine the Stent for any tumor in-growth/over-growth into the Stent lumen or whether the Stent is occluded. If the Stent lumen is clear, carefully remove using a forcep and/or snare. Grasp the retrieval string and/or collapse the proximal end of the Stent then carefully retrieve the Stent. If the Stent cannot be easily withdrawn, do not remove the Stent.

Caution: Do not allow excessive force to remove the stent as it may cause disconnection to the retrieval string.

To reposition a HILZO® Covered Stent immediately after deployment, use forceps or a snare to grasp the retrieval string and gently adjust to the correct placement.

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Storage: Store at room temperature.

Disposal Requirements:

After use, dispose product and packaging in accordance with hospital, administrative and/or local government policy for biohazardous medical waste.

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HILZO® Colonic Uncovered Stent
HILZO® Colonic Covered Stent

Figure 1 illustrates the experimental setup. A participant is seated at a table, interacting with a 'Control system' (a computer monitor and keyboard) and a 'Control panel' (a set of buttons and a joystick). The participant is also connected to a 'Control system' (a computer monitor and keyboard) and a 'Control panel' (a set of buttons and a joystick). The participant is also connected to a 'Control system' (a computer monitor and keyboard) and a 'Control panel' (a set of buttons and a joystick).

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