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IFU Riptide Aspiration System
Description of documents

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State of California


County of Orange

Subscribed and sworn to (or affirmed) before me on this 4th day of February,
2021 by Michelle Mirzaee, proved
to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

Steve Oberlander

Signature

(Seal)

Medtronic

Riptide™

Aspiration System

fr	Système d'aspiration
de	Aspirationssystem
it	Sistema per aspirazione
es	Sistema de aspiración
sv	aspirationssystem
nl	Aspiratiesysteem
pt-pt	Sistema de aspiração
fi	Aspiraatiojärjestelmä
da	Aspirationssystem
el	Σύστημα αναρρόφησης
cs	Aspirační systém
hu	Aspirációs rendszer
ru	Аспирационная система
pl	System do aspiracji
tr	Aspirasyon Sistemi
no	Aspirasjonssystem
sk	Aspiračný systém
ro	Sistem de aspirație
bg	Система за аспирация
et	Aspiratsioonisüsteem
hr	Sustav za aspiraciju
lt	aspiracinė sistema
lv	aspirācijas sistēma
sl	Aspiracijski sistem
pt-br	Sistema de aspiração
ko	흡인 시스템
ar	نظام شفط

Mulle Ball
03 FEB 2021

TABLE OF CONTENTS

Riptide™ Aspiration System

English.....	3
Français.....	4
Deutsch.....	5
Italiano.....	6
Español.....	8
Svenska.....	9
Nederlands.....	10
Português (Portugal).....	12
Suomi.....	13
Dansk.....	14
Ελληνικά.....	15
Česky.....	17
Magyar.....	18
На русском языке.....	19
Polski.....	21
Türkçe.....	22
Norsk.....	23
Slovenčina.....	25
Română.....	26
Български.....	27
Eesti.....	29
hrvatski.....	30
Lietuvių.....	31
Latviski.....	32
Slovenščina.....	33
Português (Brasil).....	34
한국어.....	36
العربية.....	38
Symbol Glossary.....	39

English Instructions for Use

Riptide™ Aspiration System

DEVICE DESCRIPTION

The Riptide™ Aspiration System is comprised of the following devices:

• React™ 68 Catheter	REACT-68
• React™ 71 Catheter	REACT-71
• Riptide™ Large Bore Aspiration Tubing	MAT-110-110
• Riptide™ Aspiration Pump	MAP-1000EU
• Riptide™ Collection Canister with Intermediate Tubing	MAC-1200

The Riptide™ Aspiration System is designed to remove occlusive thrombus from the cerebral vasculature using continuous aspiration. The catheter shaft has a hydrophilic coating that spans the distal 40 cm to reduce friction during use. The catheter acts as a conduit for providing direct aspiration to the occlusion site generated by the Riptide™ Aspiration Pump. The catheter is introduced through a guide catheter or sheath and into the intracranial vasculature and guided over a neurovascular guidewire and/or microcatheter to the site of the primary occlusion. The catheter is connected to the Riptide™ Aspiration Pump through the aspiration tubing. The Intermediate Tubing connects the Riptide™ Collection Canister to the Riptide™ Aspiration Pump. The Riptide™ Aspiration System in this configuration may be used to aspirate thrombus from the occluded vessel.

Note: The catheter and aspiration tubing are Applied Parts.

INTENDED PURPOSE / INDICATIONS FOR USE

The Riptide™ Aspiration System is intended for use in the revascularization of patients with acute ischemic stroke secondary to intracranial large vessel occlusive disease (within the internal carotid, middle cerebral M1 and M2 segments, basilar, and vertebral arteries) within 8 hours of symptom onset. Patients who are ineligible for intravenous tissue plasminogen activator (IV t-PA) or who fail IV t-PA therapy are candidates for treatment.

CONTRAINDICATION

There are no known contraindications.

WARNINGS

- The Riptide™ Aspiration System should only be used by physicians who have received appropriate training in interventional neuro-endovascular techniques for treatment of acute ischemic stroke. The Riptide™ Aspiration System is designed for setup by clinical support staff under the supervision of a trained physician.
- Do not advance, retrieve, or use any component of the Riptide™ Aspiration System against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device or system as a unit. Unrestrained torqueing or forced insertion of the catheter against resistance may result in damage to the device or vessel.
- For vessel safety, do not perform more than three attempts at aspiration in the same vessel using the Riptide™ Aspiration System.
- Do not use devices in conjunction with or in place of those devices specified above when using the Riptide™ Aspiration System.

PRECAUTIONS

- The catheter and the aspiration tubing are supplied sterile and are intended for single use only. Do not resterilize or re-use. Resterilization and/or re-use may result in ineffective catheter coating lubrication, which may result in high friction and the inability to access the target neurovascular location.
- Use the catheter, aspiration tubing, and Riptide™ Collection Canister with Intermediate Tubing prior to the "Use-by date".
- The Riptide™ Aspiration Pump and Riptide™ Collection Canister with Intermediate Tubing are supplied non-sterile and are intended for use outside of the sterile field only.
- The Riptide™ Collection Canister with Intermediate Tubing is intended for single use only. Dispose in accordance with standard biological waste disposal procedures.
- Do not use kinked or damaged devices. Do not use open or damaged packages.
- Inspect the catheter before use to verify that its size and condition are suitable for the specific procedure.
- Use the catheter under fluoroscopic guidance to confirm placement of the catheter proximal to the thrombus prior to applying aspiration.
- Maintain a constant infusion of appropriate flush solution.
- When performing aspiration, ensure that the aspiration tubing flow switch is in the open position for only the minimum time needed to remove thrombus.
- The recommended continuous duration of aspiration is a maximum of 90 seconds.
- If repositioning of the catheter is necessary during the revascularization procedure, such repositioning should be performed over an appropriate neurovascular guidewire and/or microcatheter using standard techniques.
- Do not use automated high-pressure contrast injection equipment with the catheter because it may damage the device.
- Medical management and post stroke care should follow AHA/ASA guidelines and institutional protocols.
- As in all surgical interventions, monitoring of intra-procedural blood loss is recommended so that appropriate management may be instituted.
- Possible malfunctions of the Riptide™ Aspiration System include, but are not limited to, issues associated with vacuum for the aspiration tubing, Riptide™ Aspiration Pump,

and Riptide™ Collection Canister with Intermediate Tubing that may lead to decreased therapeutic response during mechanical thrombectomy.

Operators should take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

POTENTIAL ADVERSE EVENTS

- acute occlusion
- allergic reaction and anaphylaxis from contrast media
- arrhythmia
- complication at puncture site
- complication of radiation exposure (e.g., alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia)
- death
- decreased therapeutic response
- embolism
- fistula
- foreign body in patient
- hypertension/hypotension
- infection
- inflammation
- intracranial hemorrhage
- ischemia
- nerve damage
- neurological deficits including stroke
- organ failure/shock
- pain
- stenosis
- vessel spasm, thrombosis, dissection, or perforation

PROCEDURE

1. Refer to Warnings, Precautions, and Potential Adverse Events prior to use.
2. As each component of the Riptide™ Aspiration System is used, remove the component from the packaging and inspect for damage or for kinks.
3. Prepare a guide catheter or sheath according to the device manufacturer's instructions for use.
4. Place the guide catheter or sheath into the appropriate artery for access to the occlusion site.

CATHETER PREPARATION AND USE

1. Confirm vessel diameter and ensure the catheter is appropriate for treatment of the vessel.

Recommended Sizing Guidelines	
Model	Recommended Vessel Diameter (mm)
	Minimum
REACT-68	≥2.7
REACT-71	≥2.7

2. Prepare the catheter for use by flushing the packaging hoop and device with heparinized saline.
3. Remove the catheter from the packaging hoop.
4. Inspect the catheter for kinks or damage. If any damage is observed, replace with a new device.
- NOTE:** When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).
5. Connect a rotating hemostasis valve to the proximal hub of a guide catheter or sheath.
6. Insert the catheter through the rotating hemostasis valve. If required, the provided split y-introducer may be used to facilitate the introduction of the catheter through the rotating hemostasis valve, inserting the tip through the introducer starting at the split end and penetrating the valve with the tubular introducer end. After the catheter has passed through the rotating hemostasis valve, the split y-introducer may be removed by backing it out of the valve, pulling gently on the split ends to separate the remaining portion. Discard the separated introducer.
7. Using conventional catheterization techniques under fluoroscopic guidance advance the catheter into the target vessel over an appropriate neurovascular guidewire and/or microcatheter. Position the catheter proximal to the thrombus. Remove the guidewire and/or microcatheter from the catheter.

NOTE: Access to the site of occlusion may be aided by the triaxial use of a compatible inner microcatheter.


Figure 1: Riptide™ Aspiration System Setup Diagram

Riptide™ Aspiration System PREPARATION AND USE

1. Refer to Figure 1 (above) for the correct arrangement of devices for proper setup of the Riptide™ Aspiration System.
2. Refer to the Riptide™ Aspiration Pump User's Manual for pump setup, operating instructions and troubleshooting if required.
3. Refer to the Riptide™ Collection Canister Quick Reference Guide for assembly setup.

4. Place the Riptide™ Aspiration Pump on a stable, flat horizontal surface outside of the sterile field within reach of the aspiration tubing and providing sufficient slack to the patient table. Locate the Riptide™ Aspiration Pump such that a minimum of 6 inches (approximately 15 cm) clearance exists in all directions to adjacent objects or other common cath lab equipment, powered or un-powered.
5. Attach the Riptide™ Aspiration Pump's power cord to the designated mains power specified in the Riptide™ Aspiration Pump User's Manual. Confirm that the power cord is securely plugged into the back of the pump enclosure. If required, attach the potential equalization connector to the potential equalization plug.
6. The Riptide™ Collection Canister lid should be pre-assembled to the cup as it is removed from its packaging. If it is not or becomes separated, secure the lid to the cup by snapping it tightly on the lip.
7. The large drain port on the Riptide™ Collection Canister lid should be pre-assembled with a closure cap as it is removed from its packaging. If it is not or becomes separated, seal the port with the large tethered cap on the lid by pressing it in place firmly. The small caps attached to the lid may remain tethered and unused until needed.
8. Place the Riptide™ Collection Canister into the receptacle on the Riptide™ Aspiration Pump and press firmly into place. The horizontal port marked with the embossed label as labeled "PATIENT" should point away from the pump enclosure to enable unobstructed attachment of the aspiration tubing.
9. The Intermediate Tubing should be pre-assembled to the Riptide™ Collection Canister lid as it is removed from its packaging. If it is not or becomes separated, secure the connector at one end of the supplied Intermediate Tubing over the vertical central port of the Riptide™ Collection Canister marked with the embossed label as "VACUUM" by pressing down firmly. Either connector end of the Intermediate Tubing may be used.
10. Attach other end of the Intermediate Tubing to the vacuum inlet of the Riptide™ Aspiration Pump by pressing down firmly.
11. Open the aspiration tubing pouch with the contents accessible toward the sterile field. Remove the tubing assembly from the pouch into the sterile field direction. Place or hold the distal end with the flow control switch on the patient table and pass the proximal end of aspiration tubing from the sterile field to the Riptide™ Collection Canister located on the Riptide™ Aspiration Pump.
12. Attach the connector at the proximal end of the aspiration tubing to the horizontal Riptide™ Collection Canister port marked with the embossed label as labeled "PATIENT". Press the connector firmly onto the port.
13. Slide the actuator of the flow control switch on the aspiration tubing to OFF to close the tubing to flow.
14. Press the start button on the Riptide™ Aspiration Pump and allow it to run for at least one minute. Confirm normal operation and the presence of vacuum as indicated on the gauge.
15. Using the vacuum control knob on the Riptide™ Aspiration Pump, adjust the vacuum to a reading of minimum 20 inHg (68 kPa) on the gauge but not exceeding 25 inHg (85 kPa). Adjust the vacuum by rotating the vacuum control knob until the gauge reads the desired vacuum level. To increase the vacuum, turn the knob clockwise. To decrease the vacuum, turn the knob counter-clockwise.
- NOTE:** It is normal for the indicator needle of the gauge to oscillate or vibrate in response to pump pulsation. Allow the range of the needle to stabilize and adjust the vacuum control so that the midpoint of the motion aligns with the desired setting.
16. The Riptide™ Aspiration System is now ready to aspirate.
17. Connect luer fitting at the distal end of the aspiration tubing to the side port of the rotating hemostasis valve attached to the catheter.
18. Advance the distal tip of the catheter to engage with the proximal end of the thrombus or in accordance with the physicians preferred technique.
19. To begin aspiration, slide the flow control switch toward the ON position. To stop the aspiration flow, slide the actuator of the flow control switch toward the OFF position.
20. If necessary, additional passes with catheter may be used to further remove thrombus at the discretion of the physician.
21. After withdrawal from the patient, disconnect the aspiration tubing from the catheter; replace it with 5cc or 10cc syringe, and aspirate approximately 5cc of blood from the catheter to remove any thrombus that may remain.
22. Obtain a post-treatment angiogram in accordance with the physicians preferred technique.
23. At the conclusion of the aspiration procedure, turn OFF the Riptide™ Aspiration Pump. The catheter, the aspiration tubing, and the Riptide™ Collection Canister with Intermediate Tubing are single use. Remove and discard in accordance with standard biological waste disposal procedures.
24. Clean and prepare the Riptide™ Aspiration Pump for storage as described in the User's Manual.

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COPY CERTIFICATION BY DOCUMENT CUSTODIAN

I, Michelle Mirzaee, of legal age, being first duly sworn, hereby
Name of document custodian

Swear (or affirm) that the attached reproduction of

Riptide Aspiration Pump User Manual
Description of documents

Is a true and exact copy of the correct and complete original document.

Michelle Mirzaee

Signature of document custodian (affiant)

A notary public or other officer completing this certificate verifies only the identity of the individual who signed the document to which this certificate is attached, and not the truthfulness, accuracy, or validity of that document.

State of California


County of Orange

Subscribed and sworn to (or affirmed) before me on this 18th day of November,
2021 by Michelle Mirzaee —, proved

to me on the basis of satisfactory evidence to be the person~~s~~ who appeared before me.

Steve Oberlander
Signature

(Seal)

Medtronic

Riptide™

Aspiration Pump

fr	Pompe d'aspiration
de	Aspirationspumpe
it	Pompa per aspirazione
es	Bomba de aspiración
sv	Aspirationspump
nl	Aspiratiepomp
pt-pt	Bomba de aspiração
fi	Aspiraatiopumppu
da	Aspirationspumpe
el	Αντλία αναρρόφησης
cs	Aspirační pumpa
hu	Aspirációs pumpa
ru	Аспирационная помпа
pl	Pompa do aspiracji
tr	Aspirasyon Pompası
no	Aspirasjonspumpe
sk	Aspiračná pumpa
ro	Pompă de aspirație
bg	Помпа за аспирация
et	Aspiratsioonipump
hr	Aspiracijska pumpa
lt	Siurbimo pompa
lv	aspirācijas sūkņi
sl	Aspiracijska črpalka
pt-br	Bomba de aspiração
ko	흡인 펌프
ar	مضخة الشفط

Mull Hall

18 NOV 2021

TABLE OF CONTENTS

Riptide™ Aspiration Pump

English	3
Français	7
Deutsch	11
Italiano	15
Español	19
Svenska	23
Nederlands	27
Português (Portugal)	31
Suomi	35
Dansk	39
Ελληνικά	43
Česky	48
Magyar	52
На русском языке	56
Polski	61
Türkçe	66
Norsk	70
Slovenčina	74
Română	78
Български	82
Eesti	87
Hrvatski	91
Lietuvių	95
Latviski	99
Slovenščina	103
Português (Brasil)	107
한국어	111
العربية	118
Symbol Glossary	119

English User's Manual

en

Riptide™ Aspiration Pump 220-240VAC, 50/60HZ



Figure 1: Riptide™ Aspiration Pump Features and Controls


Figure 2: Riptide™ Aspiration Pump and Riptide™ Collection Canister with Intermediate Tubing

DESCRIPTION

- The Riptide™ Aspiration Pump is an externally powered electromechanical device capable of generating a vacuum designed for drawing fluids and small particles into the Riptide™ Collection Canister during neurovascular interventional procedures. It is intended for clinical and catheter lab use and is not intended for transport or field applications.
- The recommended operating range for the Riptide™ Aspiration Pump is between 20-25 inHg (68-85 kPa).
- The Riptide™ Aspiration Pump is intended for use as a component of the Riptide™ Aspiration System. The Riptide™ Aspiration System includes the React™ 68 Catheter, React™ 71 Catheter, Riptide™ Large Bore Aspiration Tubing, Riptide™ Aspiration Pump, and Riptide™ Collection Canister with Intermediate Tubing.
- Federal law (USA) restricts this device and other components of the Riptide™ Aspiration System to sale by or on the order of a physician.
- The Riptide™ Aspiration Pump package includes the following:
 - Riptide™ Aspiration Pump
 - Region specific power cord for connection to earthed receptacle. Select correct cord from included kit.
 - Riptide™ Aspiration System Instructions for Use
 - Riptide™ Aspiration Pump User's Manual
- Components used with the Riptide™ Aspiration System (packaged separately):
 - React™ 68 Catheter, Applied Part
 - React™ 71 Catheter, Applied Part
 - Riptide™ Large Bore Aspiration Tubing, Applied Part
 - Riptide™ Aspiration Pump
 - Riptide™ Collection Canister with Intermediate Tubing

INTENDED PURPOSE / INDICATIONS FOR USE AND CONTRAINDICATIONS

- Refer to the Riptide™ Aspiration System Instructions for Use for intended purpose / indications for use and contraindications.

OPERATING DESCRIPTION


Figure 3: Riptide™ Aspiration Pump Internal Composition

The Riptide™ Aspiration Pump uses an industry proven oilless rocking piston type vacuum pump. A single, permanent split capacitor type electric motor powers a two-stage pump arrangement that provides high vacuum and efficient flow through the Riptide™ Collection Canister. In a reciprocating motion, a flexible cup mounted on top of a wristless piston and connecting rod in each stage creates vacuum as a result of the rocking motion. The vacuum is drawn in through the vacuum inlet and exhausted internally in the enclosure where it is mixed with cooling air and discharged through a vent located on the bottom.

WARNINGS

- The Riptide™ Aspiration System in contact with patients should only be used by physicians who have received appropriate training in interventional neuro-endovascular techniques and treatment of acute ischemic stroke. The Aspiration Tubing, Riptide™ Aspiration Pump, and Riptide™ Collection Canister with Intermediate Tubing are designed for setup by clinical support staff under the supervision of a trained physician.
- Do not use the Riptide™ Aspiration System with components other than the React™ 68 Catheter, React™ 71 Catheter, Riptide™ Large Bore Aspiration Tubing, Riptide™ Aspiration Pump, and Riptide™ Collection Canister with Intermediate Tubing.
- The Riptide™ Aspiration Pump is a Class I Medical Equipment. To avoid the risk of electric shock, this equipment must only be connected to a supply mains with protective earth.

PRECAUTIONS

- The Riptide™ Aspiration Pump and the Riptide™ Collection Canister with Intermediate Tubing are supplied non-sterile and are intended for use outside of the sterile field only.
- Do not use the Riptide™ Aspiration Pump without the Riptide™ Collection Canister with Intermediate Tubing connected to the vacuum inlet.
- The Riptide™ Aspiration Pump is intended for use only with the Riptide™ Collection Canister with Intermediate Tubing. The lid is equipped with integral overflow protection and a biological filter required for safe and correct operation in the Riptide™ Aspiration System. Using non-validated canisters, filters, or systems may result in improper patient connections, fluid overflow, biological contamination, and non-repairable damage of the Riptide™ Aspiration Pump.
- Do not allow aspirated fluids, overflow, or other liquids into the Riptide™ Aspiration Pump. This will cause non-repairable damage to the pump and create biological contamination. Pumps so contaminated must be removed from service and disposed of in accordance with standard biological waste disposal procedures in effect.
- Aspiration Tubing and Riptide™ Collection Canister with Intermediate Tubing are single use only.
- The Riptide™ Aspiration Pump does not contain any user serviceable parts inside. Do not attempt to open the enclosure.
- Do not use in presence of explosive atmospheres, flammable liquids or anesthetic mixtures, and nitrous oxide.
- Do not use in oxygen enriched environment.
- Do not use the Riptide™ Aspiration Pump in an environment with MRI or HF Surgical equipment.

SETUP

- Remove the Riptide™ Aspiration Pump from storage and/or its packaging (if present) and visually inspect for damage.
- Place the Riptide™ Aspiration Pump on a rigid, stable, flat, horizontal surface outside of the sterile field, within reach of the Aspiration Tubing and providing sufficient tubing slack to the patient table.
- The Riptide™ Aspiration Pump should be placed so that the tubing connections, controls and gauge are visible and accessible to the operator. The Riptide™ Aspiration Pump must never be positioned in a manner that makes it difficult to operate.
- Locate the Riptide™ Aspiration Pump such that a minimum of 6 inches (approximately 15 cm) clearance exists in all directions to adjacent objects or other common cath lab equipment, powered or un-powered. Do not stack the Riptide™ Aspiration Pump with other equipment.
- Do not block the air vents located on the back and bottom of the enclosure. Do not place the pump on thick or plush absorbent table covering and do not cover the pump with poly or other sheeting during operation. Blocked air flow will result in pump overheating.
- Attach the correct end of the supplied power cord to the power connection on the back of the enclosure. Press in firmly and confirm its solid positioning.
- If required, attach the potential equalization connector to the potential equalization plug identified by a green/yellow disk located at the back of the pump.
- Confirm that the mains power is available and the pump's nameplate ratings are compatible.

- The supplied power cords are the isolation mean and are detachable from the Riptide™ Aspiration Pump. The appliance connector is different from the mains connector. Push the cord's C13 receptacle into the appliance power inlet plug located at the back of the enclosure. At the other cord end, insert the power plug of the cord into the designated mains power receptacle – nominal 220-240VAC, 50/60HZ, 200VA minimum capacity.

NOTE: The Riptide™ Aspiration Pump should not be positioned such that it is difficult to disconnect the supplied power cord.

NOTE: Do not use electrical plug adapters instead of supplied power cord specific to user's region as use of adapter may result in damage to Riptide™ Aspiration Pump.

- The Riptide™ Aspiration Pump's power cable should not be bundled or held in proximity with other power cables that may carry very large surges and transients such as those conveyed by welders and large electric motors. In addition, separation from data cables and communications lines should follow EN50174-2:2008. In general, approximate separation from unshielded lines should be no less than 2 feet (approximately 61 cm), separation for shielded cables 2 inches (approximately 5 cm). Failure to follow these guidelines may result in electromagnetic interference due to radiated or re-radiated electrical noise.

- Check the basic operation of the Riptide™ Aspiration Pump by pressing the ON/OFF button on the front panel. The button should light up green and the pump motor should start. The presence of vacuum can be confirmed by temporarily blocking the vacuum inlet with a finger and operating the pump's vacuum control knob while observing the vacuum gauge. Note that the tethered cap is not a vacuum seal and may not seal air tight without additional finger pressure on top. Turn the pump off by pressing the ON/OFF button before proceeding.

- Do not use the Riptide™ Aspiration Pump for aspiration without the Riptide™ Collection Canister with Intermediate Tubing connected to the vacuum inlet.
- Obtain a new Riptide™ Collection Canister with Intermediate Tubing and perform the setup in accordance with the Riptide™ Aspiration System Instructions for Use and the Riptide™ Collection Canister with Intermediate Tubing Quick Reference Guide.
- The Riptide™ Aspiration Pump is now ready to connect to Aspiration Tubing per the Riptide™ Aspiration System Instructions for Use.

OPERATION

- Refer to the Riptide™ Aspiration System Instructions for Use and the Riptide™ Collection Canister with Intermediate Tubing Quick Reference Guide to confirm complete setup for aspiration.
- With the React™ 68 Catheter, or React™ 71 Catheter, Riptide™ Large Bore Aspiration Tubing, Riptide™ Aspiration Pump, and Riptide™ Collection Canister with Intermediate Tubing connected, the Riptide™ Aspiration Pump is ready to begin supplying vacuum to the connected components.
- Confirm that the flow control switch on the Aspiration Tubing is in the OFF position.
- The ON/OFF button on the control panel is a push button that requires moderate force inward to activate and de-activate. When depressed to the ON position, the button will lock in the recessed position, the green symbol will light, and the pump will begin to operate. When the button is depressed again and then released, the button will return to the original flush position, the green symbol will go dark, and pump will cease to run.
- At any time when the ON/OFF button is turned OFF after the Riptide™ Aspiration Pump is run, any vacuum present in the pump will be released. Then, the Riptide™ Collection Canister with Intermediate Tubing, and Aspiration Tubing proximal to the flow control switch will return to ambient conditions within a few seconds.
- To begin operation, press the ON/OFF button on the Riptide™ Aspiration Pump control panel and confirm that the button lights, the pump motor starts and produces a steady hum.
- Allow it to run for at least one minute. Confirm normal operation and the presence of vacuum as indicated on the gauge.
- Using the vacuum control knob on the Riptide™ Aspiration Pump, adjust the vacuum to a reading of minimum 20 inHg (68 kPa) but not exceeding 25 inHg (85 kPa) on the gauge. Adjust the vacuum by rotating the vacuum control knob. To increase the vacuum, turn the knob clockwise. To decrease the vacuum, turn the knob counter-clockwise (anti-clockwise).
- Note that it is normal for the indicator needle of the gauge to oscillate or vibrate in response to pump pulsation. Allow the range of the needle to stabilize and adjust the vacuum control knob so that the midpoint of the motion aligns with the desired setting.
- The Riptide™ Aspiration System is now ready to aspirate.
- During the procedure, monitor the vacuum level indicated on the gauge. It is normal for the needle to drift about the desired set point as use conditions change. The vacuum level may drop as a result of turning the Aspiration Tubing flow control switch on, opening a tubing connection or exchanging a catheter. The vacuum level will recover to its set point when the flow control switch is returned to OFF or when the tubing or catheter connection is restored.
- During the procedure, monitor the fluid level that may have collected in the canister. It is recommended that the fluid not be allowed to fill more than about 75% of the capacity. Should this occur, the physician should have the Riptide™ Collection Canister with Intermediate Tubing changed out.
- At the conclusion of the aspiration procedure, turn off the Riptide™ Aspiration Pump by pressing the ON/OFF button. The green light on the button will go out and the pump will stop producing a humming noise.
- The React™ 68 Catheter, React™ 71 Catheter, Riptide™ Large Bore Aspiration Tubing, and Riptide™ Collection Canister with Intermediate Tubing are single use only. Remove and discard in accordance with standard biological waste disposal procedures in effect.
- Disconnect the power cord from the mains.

CLEANING

- Place the attached tethered cap over the vacuum inlet prior to cleaning.
- Disconnect the power cord from the pump and clean separately if required.
- The outside of the Riptide™ Aspiration Pump enclosure may be cleaned by wiping down with 70% IPA and common non-bleach, non-solvent, clinical surface cleaning agents, for example Diversy Virex Tb (ready-to-use) or equivalent. Wipe dry before storage.
- Do not use cleaning agents containing petroleum base solvents, acids, caustics, or chlorinated solvents.
- The Riptide™ Aspiration Pump is not compatible with any sterilization procedures including EtO, e-beam, gamma, and autoclave.

- Do not attempt to disassemble the pump other than removing the cord and Riptide™ Collection Canister with Intermediate Tubing.
- Do not immerse the Riptide™ Aspiration Pump in water or cleaning solutions, or use a flooding rinse as a cleaning method.
- Do not inject, spray or otherwise introduce water, IPA, or cleaning agents into the cooling air intake louvers or exhaust vents located on the bottom.
- Never inject, spray, aspirate, or otherwise introduce cleaning or decontamination agents into the vacuum inlet. The pump is designed only for air contact and the introduction of any liquids or solids will cause non-repairable damage to internal components.
- If the Riptide™ Aspiration Pump is known to have ingested fluids during a procedure it should be considered as contaminated and treated in accordance with standard biological waste procedures in effect. No approved cleaning methods exist for internal decontamination.

STORAGE

- Place the tethered cap over the vacuum inlet to prevent entry of dust and other contaminants.
- The Riptide™ Aspiration Pump may be stored in a clean warehouse or storeroom commonly used for clinical and lab supplies. See characteristics below for environmental limits.
- The pump does not have an applicable "Use-by date". However, after extended storage it is recommended that it be cleaned, examined for damage, and tested for basic functionality prior to its introduction in to a clinical setting.

TROUBLESHOOTING

- No user serviceable parts are inside the pump. Do not modify the pump or open the enclosure.
- Use the troubleshooting chart below to determine common causes and recommended remedies.
- Riptide™ Aspiration Pump Fuses – replace only with exact equivalent F SA H 250V, UL listed 5 x 20 mm.

PROBLEM	POSSIBLE CAUSE	ACTION
Pump has mechanical damage to enclosure, gauge, controls	Improper handling, dropping, mechanical damage	Replace pump
Power cord or a plug on the cord damaged	Improper handling, misuse, mechanical damage	Cord may be replaced with a suitable length hospital grade or equivalent single unit meeting specifications: Mains plug: As required by site region, earthing connection required Cable: EN50525/IEC60227, integral earthing conductor required Rating: 250V, 10A Appliance socket: IEC60320-C13 Do not use plug adapters or power converters not certified to meet IEC60601-1 and IEC60601-1-2. Replacement with an inequivalent cord may decrease EMC immunity.
Power connection on pump enclosure damaged	Improper handling, misuse, mechanical damage	Replace pump
Pressing ON button causes mains breaker to trip or mains fuse to blow	Internal electrical failure	Replace pump
ON button does not light when pushed, pump does not run	Mains power not connected	Plug pump into receptacle
	Cord not connected to appliance power inlet	Plug cord into enclosure's appliance power inlet plug
	Loose cord at back of unit	Press in or adjust cord in back of enclosure
	Damaged power cord	Replace power cord
	Mains power source is bad or incorrect voltage	Plug pump into correct receptacle
	Pump fuse blown	Replace fuse(s) Use only exact equivalent, F SA H250V, UL listed 5 x 20 mm.
	Mains supply breaker not on or mains fuse blown	Turn on mains supply or replace fuse
	ON button not pushed hard enough	Press again
	Internal mechanical or electrical failure	Replace pump

PROBLEM	POSSIBLE CAUSE	ACTION
ON button does not light but pump starts and appears operating normally	Faulty ON button	Complete procedure if required but replace pump at first opportunity
ON button lights when pushed but pump makes no noise	Internal electrical failure or faulty ON button	Replace pump
ON button lights and pump starts briefly but button does not remain depressed	Faulty ON button	Replace pump
ON button lights when pushed and cooling fan on bottom runs but pump does not run	Overheated motor due to operation at full load for an extended period or blocked cooling vents	Turn pump off and allow time to cool down. Remove any absorbent padding or plastic sheeting blocking pump air vents.
	Pump contaminated by aspiration overflow	Replace pump. Units contaminated with bio fluids should be disposed of in accordance with standard hospital biohazard procedures.
	Internal mechanical or electrical failure	Replace pump
Pump runs but has abnormal operation: unusually loud, high vibration, running very hot, burning smell, smoke	Pump rocking - not sitting on level surface with all 5 supporting feet	Re-position pump on flat, rigid, horizontal surface
	Cooling air intake and/or bottom exhaust is blocked	Re-position pump. Do not place pump on plush absorbent padding and do not cover with splash sheeting.
	Pump contaminated by overflow	Replace pump. Units contaminated with bio fluids should be disposed of in accordance with standard hospital biohazard procedures.
	Pump connected to improper mains voltage supply	Replace pump with unit that matches mains power input
	Internal electrical or mechanical damage or failure	Replace pump
Gauge needle vibrates excessively	Vacuum pressure set too low – especially below 15 inHg (51 kPa).	Increase vacuum between 20-25 inHg (68-85 kPa). It is normal for the indicator needle of the gauge to oscillate or vibrate in response to pump pulsation, especially at minimal vacuums. Allow the range of the needle to stabilize and adjust the vacuum control so that the midpoint of the motion aligns with the desired setting.
	Faulty gauge if needle does not stabilize at full vacuum.	Procedure may be continued at physicians discretion but replace pump at first opportunity
Gauge needle indicates vacuum greater than 30 inHg (102 kPa) or is pegged fully anti-clockwise when running	Faulty or damaged gauge.	Use vacuum control knob to attempt to set vacuum to recommended range of 20-25 inHg (68-85 kPa). Procedure may be continued at physicians discretion but replace pump at first opportunity
Pump runs but no or low vacuum indicated at gauge, vacuum may be present at Riptide™ Large Bore Aspiration Tubing	Vacuum control knob not adjusted	Adjust knob to set vacuum to desired level

PROBLEM	POSSIBLE CAUSE	ACTION
	Faulty gauge	Use vacuum control knob to attempt to set vacuum to recommended range of 20-25 inHg (68-85 kPa). Procedure may be continued at physicians discretion but replace pump at first opportunity
	Flow control switch on Riptide™ Large Bore Aspiration Tubing not set correctly	Adjust flow control switch to OFF
	Pump contaminated by overflow, pump seizing	Replace pump
	Internal electrical or mechanical damage or failure	Replace pump
Pump runs with vacuum indicated at gauge, but no vacuum at Riptide™ Large Bore Aspiration Tubing distal end	Loose, leaking, or missing Intermediate Tubing	Install tubing correctly and check fittings
	Flow control switch in OFF position	Move switch to ON position
	Riptide™ Collection Canister lid loose or drain port open	Tighten lid, close port with tethered cap
	Riptide™ Collection Canister not vertical and secured in Riptide™ Aspiration Pump receptacle	Turn off the Riptide™ Aspiration Pump and place the Riptide™ Collection Canister in the receptacle, and restart the Riptide™ Aspiration Pump
	Riptide™ Collection Canister bio-filter plugged	Replace the Riptide™ Collection Canister. Do not attempt to clean or bypass filter.
	Plugged or kinked tubing	Straighten or replace tubing
	Broken canister, defective connector or vacuum tubing	Replace with a new Riptide™ Collection Canister with Intermediate Tubing and Riptide™ Large Bore Aspiration Tubing
	Riptide™ Collection Canister full or foaming. Overflow protection valve has activated and stopped flow	Replace the Riptide™ Collection Canister. Do not attempt to bypass overflow valve.
	Riptide™ Large Bore Aspiration Tubing flow control switch set incorrectly	Set flow switch to desired setting
	Defective or broken Riptide™ Large Bore Aspiration Tubing flow control switch	Replace the Riptide™ Large Bore Aspiration Tubing
Pump runs but low suction or insufficient vacuum in Catheter	Loose, leaking, or missing Intermediate Tubing	Install tubing and check fittings
	Loose canister lid, vacuum leak	Reposition and press lid closed
	Canister dump/vent port is open or leaking	Close dump/vent port with tethered cap on lid
	Pinched or kinked connector or suction tubing	Straighten connector or suction tubing
	Incorrect vacuum control knob setting	Adjust vacuum control knob to proper setting
	Incorrect flow control switch setting	Adjust flow control switch to proper setting – OFF to build vacuum, ON to direct vacuum to catheter
	Flow control switch on tubing defective	Replace the Riptide™ Large Bore Aspiration Tubing
	Canister bio filter is clogged or overflow valve is closed	Replace the Riptide™ Collection Canister. Do not attempt to clean or bypass filter.
	Vacuum leak internal to pump enclosure	Replace pump
	Canister is full and overflow protection working normally	Replace the Riptide™ Collection Canister. Do not attempt to clean or bypass canister. Dispose of full canister in accordance with standard hospital biohazard procedures.

PROBLEM	POSSIBLE CAUSE	ACTION
	Pump contaminated by canister overflow or liquid ingress from previous procedure	Replace pump. Do not attempt to clean or decontaminate pump. Pump units contaminated with bio fluids should be disposed of in accordance with standard hospital biohazard procedures.
Vacuum indicated on gauge does not respond to change in vacuum control setting	Faulty vacuum control knob or faulty vacuum gauge	Replace pump
	Internal electrical or mechanical damage or failure	Replace pump
Pump spontaneously stops running or is intermittent during use	Loss of mains power, breaker interruption, power surge, brownout	Restore mains power or connect to uninterruptible power supply.
	Power cord disturbed, loose	Re-plug power cord securely to mains and pump enclosure
	Blown Riptide™ Aspiration Pump fuse	Replace fuse(s) Use only exact equivalent F 5A H 250V, UL listed 5 x 20 mm.
	Pump overheated due to long continuous running in poor cooling	Allow pump to cool, correct cooling air flow blockage, and restart
	Riptide™ Collection Canister overflowed and fluid ingress has seized pump	Replace Riptide™ Aspiration Pump and Riptide™ Collection Canister
	Attempted use of non-validated canister without filters; fluid ingress causing pump to seize	Replace Riptide™ Aspiration Pump. Use correct Riptide™ Collection Canister. Pump may be biological hazard
	Incorrect tubing connection for the Riptide™ Large Bore Aspiration Tubing connected directly to vacuum inlet of the Riptide™ Aspiration Pump bypassing the Riptide™ Collection Canister	Replace pump. Pump may be biological hazard
Pump does not shut off when button is pushed	Button not pushed hard enough	Press again
	Faulty button or internal electrical failure	Pull mains plug to un-power. Replace pump
Riptide™ Collection Canister does not fit in the receptacle of the Riptide™ Aspiration Pump	Attempted use of non-validated canister system.	Replace with Riptide™ Collection Canister with Intermediate Tubing
Tubing does not fit correctly on Riptide™ Aspiration Pump vacuum inlet	Attempted use of non-validated canister and/or tubing	Replace with Riptide™ Collection Canister with Intermediate Tubing
	Attempted connection of the Riptide™ Large Bore Aspiration Tubing directly to the vacuum inlet of the Riptide™ Aspiration Pump	Review the Riptide™ Aspiration System Instructions for Use. Install the Riptide™ Collection Canister and connect the Riptide™ Large Bore Aspiration Tubing to the port labeled "PATIENT" on the Riptide™ Collection Canister.

MAINTENANCE

- Operating life: The Riptide™ Aspiration Pump has an operating life of 500 hours.
- The enclosure does not contain any user serviceable parts inside and should not be opened.
- Fuses: Located at the back of the unit. Replace only with exact equivalent F 5A H 250V, UL listed 5 x 20mm.
- Sterilization: The Riptide™ Aspiration Pump is NOT designed to be sterilized by autoclaving, chemical means (such as ETO and other sterilant compounds), or radiation. Any attempt to do so will cause non-repairable damage to the unit.

DISPOSAL

- Riptide™ Aspiration Pump with mechanical/electrical failure or normal wear-out: Dispose of non-functioning pumps in accordance with standard hospital procedures for electromechanical equipment.
- Riptide™ Aspiration Pump with contamination by aspirated fluid ingress: Dispose of in accordance with standard hospital biohazard procedures.

* The React™ 68 Catheter, React™ 71 Catheter, Riptide™ Large Bore Aspiration Tubing and Riptide™ Collection Canister with Intermediate Tubing are single use only. Dispose of in accordance with standard hospital biohazard procedures.

CHARACTERISTICS

Characteristic	Value
Model	Riptide™ Aspiration Pump MAP-1000EU
Vacuum Range, Nominal (at sea level)	0-29 inHg (0-98 kPa)
Vacuum Gauge	+/- 5% of full scale value
Flow Rate, Nominal (at sea level)	50 Hz: 0-0.66 SCFM (0-19 LPM) 60 Hz: 0-0.80 SCFM (0-23 LPM)
Power Input Rating	220-240 VAC, 50/60 Hz, 200 VA max.
Overcurrent Protection (Fuses)	F 5A H 250V, UL listed 5 x 20 mm
Approximate Dimensions	Length – 16.1 in (40.9 cm) Depth – 13.2 in (33.5 cm) Height – 12.3 in (31.2 cm)
Weight	Approximately 23 lbs (10 kg)
Duty Cycle	Non-continuous 97% (58.2 minutes on, 1.8 minutes off)
Operating Conditions (unpowered)	Temperature: 65 to 75°F (18 to 24°C) Humidity: 20-80% non-condensing Air pressure: 24-31 inHg (81-105 kPa)
Storage Conditions	Temperature: -25 to 125°F (-32 to 52°C) Humidity: 5-90% non-condensing Air pressure: 20-31 inHg (68-105 kPa), safe in standard pressurized cargo air transport
Enclosure Ingress Rating	IP21 (Drip Tight)
Power Cords	2.5 m, Unshielded, IEC Plug Types F, G, H, I, J, K, L, M, & N

ELECTRICAL SAFETY

General	IEC60601-1:2012 General standard, Medical devices
Earthbond	IEC60601-1:2012 0.2 ohm max.
Leakage Current	IEC60601-1:2012 Type BF, Class I Medical Equipment

ELECTROMAGNETIC COMPATIBILITY

The Riptide™ Aspiration Pump is expected to be relatively unaffected by common electromagnetic (EM) disturbances. In the event of unusual or extraordinary EM events beyond the limits, it is possible that device may exhibit a loss of Essential Performance (the ability to generate the necessary vacuum required for aspiration as set by the user). Should this occur, discontinue use of the device. The Riptide™ Aspiration Pump is intended for use in electromagnetic environments as specified below.

Guidance and Manufacturer's Declarations

Guidance and Manufacturer's Declaration – Electromagnetic Emissions

The Riptide™ Aspiration Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the pump should make sure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The Riptide™ Aspiration Pump may use RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B at 230 VAC, 50 Hz Class A at 220 VAC, 60 Hz	The Riptide™ Aspiration Pump is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic Emissions IEC 61000-3-2	Complies	
Voltage Fluctuations/ Flicker Emissions IEC 61000-3-3	Complies	

Guidance and Manufacturer's Declaration – Electromagnetic Immunity

The Riptide™ Aspiration Pump is intended for use in the electromagnetic environment specified below. The customer or the user of the pump should make sure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic Discharge (ESD) IEC 61000-4-2	±8 kV contact ±2 kV, ±4 kV, ±8 kV, ±15 kV air	Complies	The Riptide™ Aspiration Pump should not be affected by electrostatic discharge that might occur under normal conditions of use. It is recommended that floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2 kV 100 kHz repetition frequency	Complies	Mains power quality should be similar to that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Complies	Mains power quality should be similar to that of a typical commercial or hospital environment.
Voltage Dips, Short Interruptions, and Voltage Variations on Power Supply Lines IEC 61000-4-11	0% U _n for 0.5 cycle, 1 cycle, and 250/300 cycles, 70% U _n for 25/30 cycles	Complies	Mains power quality should be that of a typical commercial or hospital environment. If the user of the ME EQUIPMENT requires continued operation during power mains interruptions, it is recommended that the ME EQUIPMENT be powered from an uninterruptible power supply or a battery.
Power Frequency (50/60 Hz) Magnetic Field IEC 61000-4-8	30 A/m	Complies	Mains power quality should be similar to that of a typical commercial or hospital environment.

NOTE: U_n is the a.c. mains voltage prior to application of the test level.

Conducted RF IEC 61000-4-6	3V rms 150 kHz to 80 MHz	Complies	Portable RF communications equipment should be used no closer than 30 cm (12 in) to any part of the aspiration pump, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result. Interference may occur in the vicinity of equipment marked with the following symbol:
Radiated RF IEC 61000-4-3	3V/m 80 MHz to 2.7 GHz	Complies	

NOTE 1: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption, and reflection from structures, objects and people.

Recommended separation distances between portable and mobile RF communications equipment and the Riptide™ Aspiration Pump

This section provides information on the recommended separation distance between portable and mobile RF communications equipment and the Riptide™ Aspiration Pump. The Riptide™ Aspiration Pump is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the pump users can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the pump, according to the maximum output power of the communications equipment.

WARRANTY DISCLAIMER

Although this product has been manufactured under carefully controlled conditions, the manufacturer has no control over the conditions under which this product is used. The manufacturer therefore disclaims all warranties, both expressed and implied, with respect to the product including, but not limited to, any implied warranty of merchantability or fitness for a particular purpose. The manufacturer shall not be liable to any person or entity for any medical expenses or any direct, incidental or consequential damages caused by any use, defect, failure or malfunction of the product, whether a claim for such damages is based upon warranty, contract, tort or otherwise. No person has any authority to bind the manufacturer to any representation or warranty with respect to the product. The exclusions and limitation set out above are not intended to, and should not be construed so as to contravene mandatory provisions of applicable law. If any part or term of this Disclaimer of Warranty is held illegal, unenforceable or in conflict with applicable law by a court of competent jurisdiction, the validity of the remaining portions of this Disclaimer of Warranty shall not be affected, and all rights and obligations shall be construed and enforced as if this Disclaimer of Warranty did not contain the particular part or term held to be invalid.

Recommended Separation Distances Between RF Equipment and the Riptide™ Aspiration Pump						
Test frequency (MHz)	Band ^{a)} (MHz)	Service ^{a)}	Modulation ^{b)}	Maximum Power (W)	Distance (m)	IMMUNITY TEST LEVEL (V/m)
385	380 - 390	TETRA 400	Pulse modulation ^{b)} 18 Hz	1,8	0,3	27
450	430 - 470	GMRS 460, FRS 460	FM ^{c)} ± 5 kHz deviation 1 kHz sine	2	0,3	28
710	704 - 787	LTE Band 13, 17	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
745						
780						
810	800 - 960	GMS 800/900, TETRA 800, IDEN 820, CDMA 850, LTE Band 5	Pulse modulation ^{b)} 18 Hz	2	0,3	28
870						
930						
1 720	1 700 - 1 990	GMS 1800; CDMA 1900; GSM 1900; DECT; LTE Band 1, 3, 4, 25; UMTS	Pulse modulation ^{b)} 217 Hz	2	0,3	28
1 845						
1 970						
2 450	2 400 - 2 570	Bluetooth, WLAN, 802.11 b/g/n, RFID 2450, LTE Band 7	Pulse modulation ^{b)} 217 Hz	2	0,3	28
5 240	5 100 - 5 800	WLAN, 802.11 a/n	Pulse modulation ^{b)} 217 Hz	0,2	0,3	9
5 550						
5 780						
NOTE If necessary to achieve the IMMUNITY TEST LEVEL, the distance between the transmitting antenna and the Riptide™ Aspiration System may be reduced to 1 m. The 1 m test distance is permitted by IEC 61000-4-3.						
a) For some services, only the uplink frequencies are included. b) The carrier shall be modulated using a 50 % duty cycle square wave signal. c) As an alternative to FM modulation, 50 % pulse modulation at 18 Hz may be used because while it does not represent actual modulation, it would be worst case.						

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State of California

County of Orange


Subscribed and sworn to (or affirmed) before me on this 27th day of March,

2020 by Michelle Mirzaee, proved

to me on the basis of satisfactory evidence to be the person(s) who appeared before me.

[Signature]
Signature

(Seal)

Medtronic

ReactTM Catheter

fr	Cathéter
de	Katheter
it	Catetere
es	Catéter
sv	Kateter
nl	Katheter
pt-pt	Cateter
fi	Katetri
da	Kateter
el	Καθετήρας
cs	Katétr
hu	Katéter
ru	Катетер
pl	Cewnik
tr	Kateter
no	Kateter
sk	Katéter
ro	Cateter
bg	Катетър
et	Kateeter
hr	Kateter
lt	Kateteris
lv	katetrs
sl	Kateter
pt-br	Cateter
ko	카테터
ar	قسطرة

DWGS71009-001 Rev. C (04/2019)

TABLE OF CONTENTS

React™ Catheter

English	3
Français	4
Deutsch	5
Italiano	6
Español	7
Svenska	8
Nederlands	9
Português (Portugal)	10
Suomi	11
Dansk	12
Ελληνικά	13
Česky	14
Magyar	15
На русском языке	16
Polski	17
Türkçe	18
Norsk	19
Slovenčina	20
Română	21
Български	22
Eesti	23
hrvatski	24
Lietuvių	25
Latviski	26
Slovenščina	27
Português (Brasil)	28
한국어	29
العربية	30
Symbol Glossary	31

English Instructions for Use

en

React™ Catheter

DEVICE DESCRIPTION

The catheter is a single lumen, flexible, variable stiffness composite catheter. The catheter shaft has a hydrophilic coating that spans the distal 40 cm to reduce friction during use. The catheter shaft is visible under fluoroscopy. The catheter dimensions are included on the individual device label. The catheter is introduced through a guide catheter or sheath and into the intracranial vasculature and guided over a neurovascular guidewire and/or microcatheter. The proximal end of the catheter has a luer fitting to allow attachment of accessories and infusion of liquids through the system. The catheter is provided sterile, non-sterile, and is intended for single use only.

INDICATIONS FOR USE

The React™ 68 Catheter and React™ 71 Catheter are indicated for the introduction of interventional/diagnostic devices into the peripheral and neurovasculature.

The React™ 68 Catheter and React™ 71 Catheter are also indicated for the removal/aspiration of emboli and thrombi from selected blood vessels in the arterial system, including the neurovasculature.

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COMPATIBILITY

Refer to product label for device dimensions. Refer to labeling provided with other medical technologies to determine compatibility.

WARNINGS

- The catheter should only be used by physicians who have received appropriate training in interventional techniques.
- Do not reuse. The device is intended for single use only. Discard the catheter after one procedure. Structural integrity and/or function may be impaired through reuse or cleaning. Catheters are extremely difficult to clean after exposure to biological materials and may cause adverse patient reactions if reused.
- The catheter has not been tested for use with automated high-pressure contrast injection equipment, do not use this equipment with the device because it may damage the device.

PRECAUTIONS

- Do not use open, kinked, or damaged devices.
- Do not autoclave.
- Use the catheter in conjunction with fluoroscopic visualization.
- Inspect the catheter before use to verify that its size and condition are suitable for the specific procedure.
- Do not advance or withdraw the catheter against resistance without careful assessment of the cause using fluoroscopy. If the cause cannot be determined, withdraw the device. Moving or torquing the device against resistance may result in damage to the vessel or device.
- Maintain a constant infusion of appropriate flush solution.
- If flow through the catheter becomes restricted, do not attempt to clear the lumen by infusion. Remove and replace the device.
- Extreme care must be taken to avoid damage to the vasculature through which the catheter passes. The catheter may occlude smaller vessels. Care must be taken to avoid complete blood flow blockage.
- Torquing the catheter may cause damage which could result in kinking and possible separation along the catheter shaft. Should the catheter become severely kinked, withdraw the catheter.
- An appropriate anticoagulation therapy should be applied per institutional guidelines.
- Operators should take all necessary precautions to limit X-radiation doses to patients and themselves by using sufficient shielding, reducing fluoroscopy times, and modifying X-ray technical factors where possible.

POTENTIAL ADVERSE EVENTS

Possible complications include, but are not limited to, the following:

- acute occlusion
- allergic reaction and anaphylaxis from contrast media
- complication at puncture site
- complication of radiation exposure (e.g., alopecia, burns ranging in severity from skin reddening to ulcers, cataracts, and delayed neoplasia)
- embolism
- false aneurysm formation
- infection
- inflammation
- intracranial hemorrhage
- ischemia
- neurological deficits including stroke
- vessel spasm, thrombosis, dissection or perforation

RECOMMENDED PREPARATION AND USE

1. Gently remove the catheter and packaging hoop from the pouch by grasping the packaging hoop and slowly pulling the catheter and hoop out of the pouch.
2. Flush package hoop with saline before removing catheter.
3. Remove the catheter from the packaging hoop.
4. Inspect the catheter for kinks or other damage. If any damage is observed, replace with new device.
5. Connect a rotating hemostasis valve (RHV) to the hub of the catheter and flush the inner lumen with heparinized saline.
6. Introduce the catheter into the vasculature through the sheath, and over a guide wire using a percutaneous entry technique of choice. A split sheath is provided in the packaging to provide support and facilitate the introduction of the tip of the catheter into the sheath.
NOTE: The maximum diameter recommended for a guide wire is 0.038" (0.97 mm). When using the React™ 71 Catheter, the minimum inner diameter recommended for a sheath is 0.087" (2.21 mm).
7. Under fluoroscopic guidance, advance the system over the guide wire until the desired position is attained.
8. Remove the guide wire.
9. If introducing interventional devices, refer to the instructions for use provided by the manufacturer.
10. If performing aspiration, apply aspiration to the catheter.
11. If necessary, additional attempts at flow restoration may be used.
12. Obtain a post-treatment angiogram in accordance with the physicians preferred technique.

STORAGE

Store in a cool, dry and dark place. See the product label for the device "Use-by date". Do not use the device beyond its labeled "Use-by date".

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© 2019 Medtronic.
Medtronic, Medtronic with logo and Medtronic logo are U.S. and internationally registered trademarks of Medtronic.
Therapeutics, Inc. Dba Ev3