COVIDIEN

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Covidien llc, USA

Sr. Regulatory Affairs Operations Specialist

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"Operating documents" for submission in the Russian Federation:

Electrosurgical Monopolar Electrode Valleylab in modifications.

Manufactured by:

Covidien llc, USA. 15 Hampshire Street, Mansfield MA 02048, USA

Manufacturing site:

Covidien 60 Middletown Avenue, North Haven, CT 06473, USA

Covidien LLC, 5920 Longbow Drive, Boulder Colorado, 80301, USA

Covidien, Boulevard Insurgentes, 19030 Libramiento, 22225 Tijuana, B. C., MEXICO

Venusa de Mexico S. de R.L. de C. V. a Lake Region Medical Company, Calle Hertz 1525, Pargue Industrial Antonio J. Bermudez, Chihuahua, 32470 Ciudad Juarez, Mexico.


State of Colorado)
) §
County of Boulder)

Subscribed and sworn to before me this 30th day of September 2021. Proved to me on the basis of satisfactory evidence to be the person whose name is subscribed to the within instrument, and acknowledged to me that she executed the same for purposes therein stated.

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For use with a max peak voltage of 5600 V.

Indications for Use

The Valleylab electrodes are indicated for use in surgical procedures (general, neurosurgical, laparoscopic, orthopedic, gynecologic, etc.) where monopolar electrosurgical cutting and coagulation are normally used.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection, or product failure risks to the patient.

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- Oxygen-enriched environments
- Oxidizing agents such as nitrous oxide (N₂O) atmospheres
 - Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
 - If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.
- Alcohol-based skin prepping agents and tinctures
 - Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.

Warning

- Naturally occurring flammable gases (such as methane) that may accumulate in body cavities

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- When using electrosurgery in the same room with gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes or near the surgical site.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen-enriched environments such as in throat or mouth procedures. Eschar plus high oxygen may create embers. Keep the electrode clean and free of all debris.
- Facial and other body hair is flammable. Water-soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

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Warning

Some surgeons may elect to “buzz the hemostat” during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon’s hands may result. To minimize the risk, take these precautions:

- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.
- “Buzz the hemostat” below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon’s hands.

Warning

- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.

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The needles are fragile and should be handled with care to avoid damage to the needle and injury to hospital personnel.

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Activate the electrosurgical unit **only** when ready to deliver electrosurgical current and the active tip is in view (especially if looking through an endoscope).

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The electrodes are intended for single use only. Safely **discard after use** to prevent injury to hospital personnel. These electrodes are not designed to withstand resterilization. **Do not resterilize.**

Important

Wipe the electrode often with moist gauze or other material.

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

Instructions for Use

Note: The electrodes can be used with both disposable and reusable handswitching and footswitching electrosurgical pencils.

1. Ensure the pencil is not connected to the generator.


2. Grasp the insulating sleeve on the electrode. Insert the electrode into the pencil.


3. Ensure the insulating sleeve fits securely inside the nose of the pencil so that the nose overlaps the insulating sleeve by at least 1/8" (0.3 cm).

4. Hex electrodes have a depth indicator. The line on the depth indicator should be flush with the tip of the handswitching pencil.

5. A tip protector covers the end of some electrodes. If a tip protector is present, remove it before use.

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
0086

Part No. 1076247

COVIDIEN, COVIDIEN with logo, and Covidien logo and Positive Results for Life are U.S. and internationally registered trademarks of Covidien AG. Other brands are trademarks of a Covidien company, TM brands are trademarks of their respective owner.

© 2012 Covidien.

 Covidien llc,
15 Hampshire Street, Mansfield, MA 02048 USA.

 Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Printed in USA.

www.covidien.com

REV 08/2014

STERILE EO

**Rx
ONLY**

Consult
instructions
for use

Do not reuse

MEAD

0086

Part No. 1056398

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Other brands are trademarks of a Covidien company.

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0086

Part No. 1056398

COVIDIEN, COVIDIEN with logo and Covidien logo are U.S. and internationally registered trademarks of Covidien AG.

Other brands are trademarks of a Covidien company.

© 2012 Covidien.

Covidien Inc,
15 Hampshire Street, Mansfield, MA 02048 USA.
Covidien Ireland Limited,
IDA Business & Technology Park, Tullamore.

Printed in USA.

www.covidien.com

REV 11/2012

en

Valleylab™

- REF E1650 Straight Micro-Needle Electrode
2 cm
- REF E1651 Straight Micro-Needle Electrode
3 cm
- REF E1652 Angled Micro-Needle Electrode
3 cm (45°)
- REF E1653 Angled Micro-Needle Electrode
3 cm (45°)
- REF E1654B Insulated Micro-Needle Electrode
9 cm

Do not use if package is opened or damaged.

2

en

Valleylab™

- REF E1650 Straight Micro-Needle Electrode
2 cm
- REF E1651 Straight Micro-Needle Electrode
3 cm
- REF E1652 Angled Micro-Needle Electrode
3 cm (45°)
- REF E1653 Angled Micro-Needle Electrode
3 cm (45°)
- REF E1654B Insulated Micro-Needle Electrode
9 cm

Do not use if package is opened or damaged.

For use with a max peak voltage of 2500 V.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Danger: Explosion Hazard Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard Both oxygen (O₂) and nitrous oxide (N₂O) support combustion. Watch for enriched O₂ and N₂O atmospheres near the surgical site, especially during head and neck surgery. Enriched O₂ atmospheres may result in fires and burns to patients or surgical personnel.

For use with a max peak voltage of 2500 V.

Warning

This product cannot be adequately cleaned and/or sterilized by the user in order to facilitate safe reuse, and is therefore intended for single use. Attempts to clean or sterilize these devices may result in a bio-incompatibility, infection or product failure risks to the patient.

Danger: Explosion Hazard Do not use electrosurgery in the presence of flammable anesthetics.

Fire Hazard Both oxygen (O₂) and nitrous oxide (N₂O) support combustion. Watch for enriched O₂ and N₂O atmospheres near the surgical site, especially during head and neck surgery. Enriched O₂ atmospheres may result in fires and burns to patients or surgical personnel.

Warning

Fire/Explosion Hazard The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen enriched environments
- Oxidizing agents such as nitrous oxide (N₂O) atmospheres
 - Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
 - If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.

Warning

- Alcohol-based skin prepping agents and tinctures
 - Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
- Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

Fire Hazard The sparking and heating associated with electrosurgery can provide an ignition source.

Observe fire precautions at all times:

- When using electrosurgery in the same room with gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes or near the surgical site.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments such as in throat or mouth procedures. Eschar plus high oxygen may create embers. Keep the electrode clean and free of all debris.

3

Warning

Fire/Explosion Hazard The following substances contribute to increased fire and explosion hazards in the operating room:

- Oxygen enriched environments
- Oxidizing agents such as nitrous oxide (N₂O) atmospheres
 - Verify all anesthesia circuit connections are leak free before and during use of electrosurgery.
 - Verify endotracheal tubes are leak free and that the cuff seals properly to prevent oxygen leaks.
 - If an uncuffed tube is in use, pack the throat with wet sponges around the uncuffed tube.
 - If possible, stop supplemental oxygen at least one minute before and during use of electrosurgery.

Warning

- Alcohol-based skin prepping agents and tinctures
 - Activate the electrosurgical unit only after vapors from skin prep solutions and tinctures have dissipated.
- Naturally occurring flammable gases (such as methane) that may accumulate in body cavities.

Fire Hazard The sparking and heating associated with electrosurgery can provide an ignition source.

Observe fire precautions at all times:

- When using electrosurgery in the same room with gases or flammable substances, prevent pooling of fluids and the accumulation of gases under surgical drapes or near the surgical site.
- Tissue buildup (eschar) on the tip of an active electrode poses a fire hazard, especially in oxygen enriched environments such as in throat or mouth procedures. Eschar plus high oxygen may create embers. Keep the electrode clean and free of all debris.

Warning

- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

The electrode must fit completely and securely into the pencil. An incorrectly seated electrode may result in burns to the patient or surgical personnel.

Fire Hazard Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place longer electrodes such as extended electrodes away from the patient and drapes.

Warning

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

Some surgeons may elect to "buzz the hemostat" during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon's hands may result. To minimize the risk, take these precautions:

- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.

Warning

- Facial and other body hair is flammable. Water soluble surgical lubricating jelly may be used to cover hair close to the surgical site to decrease flammability.

The electrode must fit completely and securely into the pencil. An incorrectly seated electrode may result in burns to the patient or surgical personnel.

Fire Hazard Always place the active electrode in a clean, dry, insulated safety holster when not in use.

- Electrosurgical accessories that are activated or hot from use can cause unintended burns to the patient or surgical personnel.
- Electrosurgical accessories may cause fire or burn if placed close to or in contact with flammable materials such as gauze or surgical drapes. Place longer electrodes such as extended electrodes away from the patient and drapes.

Warning

Confirm proper electrosurgical settings prior to and during a procedure. Use the lowest power settings to achieve the desired effect. If increased power settings are requested, check the patient return electrode and all accessory connections before major power setting adjustments.

Some surgeons may elect to "buzz the hemostat" during surgical procedures. It is not recommended, and the hazards of such a practice probably cannot be eliminated. Burns to the surgeon's hands may result. To minimize the risk, take these precautions:


- Do not lean on the patient, the table, or the retractors while buzzing the hemostat.
- Activate cut rather than coag. Cut has a lower voltage than coag.
- Use the lowest power setting possible for the minimum time necessary to achieve hemostasis.
- Activate the generator after the accessory makes contact with the hemostat. Do not arc to the hemostat.

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
Warning

- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.
- "Buzz the hemostat" below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.
- "Buzzing hemostats" may damage the electrode.

Caution

Do not exceed maximum power limits as stated in Instructions for Use. Exceeding these power settings may result in patient injury or product damage.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Caution

Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages.

Needle electrodes will easily puncture surgical gloves or tissue. Handle and dispose of electrodes carefully.

Needle electrodes are fragile. Handle them with care to avoid damage to the needle and injury to hospital personnel.

Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors.

Do not modify or add to the insulation of active electrodes.

5

Warning

- Firmly grasp as much of the hemostat as possible before activating the generator. This disperses the current over a larger area and minimizes the current concentration at the fingertips.
- "Buzz the hemostat" below hand level (as close as possible to the patient) to reduce the opportunity for current to follow alternate paths through the surgeon's hands.
- "Buzzing hemostats" may damage the electrode.

Caution

Do not exceed maximum power limits as stated in Instructions for Use. Exceeding these power settings may result in patient injury or product damage.

Always use the lowest power setting that achieves the desired surgical effect. Use the active electrode for the minimum time necessary in order to reduce the possibility of unintended burn injury.

Caution

Pediatric applications and/or procedures performed on small anatomic structures may require reduced power settings. The higher the current flow and the longer the current is applied, the greater the possibility of unintended thermal damage to tissue, especially during use on small appendages.

Needle electrodes will easily puncture surgical gloves or tissue. Handle and dispose of electrodes carefully.

Needle electrodes are fragile. Handle them with care to avoid damage to the needle and injury to hospital personnel.

Before use, examine the electrosurgical unit and accessories for defects. Do not use cables or accessories with damaged (cracked, burned, or taped) insulation or connectors.

Do not modify or add to the insulation of active electrodes.

Caution

Activate the electrosurgical unit **only** when you are ready to deliver electrosurgical current and the active tip is in view (especially if looking through an endoscope).

Deactivate the electrosurgical unit **before** the tip leaves the surgical site.

The electrodes are intended for single use only. Safely **discard after use** to prevent injury to hospital personnel. These electrodes are not designed to withstand resterilization. Do not resterilize.

Notice

Needle electrodes are designed for precise low power use during monopolar electrosurgery. Using a needle at high power settings for extended periods of time may result in damage to the needle. Use low power settings for short periods of time to prevent needle damage.

Do not bend the electrode. Product damage may result.

Notice

Electrosurgical generators (e.g., Force FX™ or Force EZ™ generators) produce desired surgical effects at lower power cut mode settings than conventional electrosurgical generators.

Important

Keep the electrode clean. Wipe the electrode often with moist gauze or similar material.

Instructions for Use**Notice**

The maximum recommended power for micro-needle electrodes is **25 watts**.

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

3. Ensure the insulating sleeve fits securely inside the nose of the pencil so that the nose overlaps the insulating sleeve by at least 1/8" (0.3 cm).
4. For the E1654B insulated micro-needle electrode, make sure the insulating sleeve fits securely over the nose to the indicator line.

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
3. Ensure the insulating sleeve fits securely inside the nose of the pencil so that the nose overlaps the insulating sleeve by at least 1/8" (0.3 cm).
4. For the E1654B insulated micro-needle electrode, make sure the insulating sleeve fits securely over the nose to the indicator line.

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