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
1. That the attached document(s) Appendix to the Instructions For Use
SALUS-RSWT-DUET

☒ is(are) true and correct.

☐ is(are) true translation for the text(s) originally written in the Korean
language conscientiously believing the same to be true and correct.

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Appendix to the Operating Manual

Shock wave therapy (SWT) device, model SALUS-RSWT-DUET

(hereinafter referred to as UVT device(s), product(s), device(s), SALUS-RSWT-DUET)

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SAFETY INSTRUCTIONS

ELECTRICAL SAFETY PRECAUTIONS:

- The rated power supply voltage for the UVT device is 230 V, alternating current;
- Make sure that all connecting parts (power cord or component cables) are correctly connected to the UVT device;
- Make sure that the UVT device is completely grounded;
- Repairs, extensions and installations of the UVT apparatus must not be carried out by anyone other than specialized personnel authorized by the manufacturer. Arbitrary disassembly/assembly of the UVT device by the user is strictly prohibited;
- Before connecting other equipment not specified in the instruction manual, be sure to consult the manufacturer or the manufacturer's authorized representative;
- To avoid electrical interference during use, the UVT device should be installed at a considerable distance from the power source, any generator, X-ray equipment, transmitting device, mobile electrical wire, etc.;
- Be careful because of possible electromagnetic or other interference between the UVT device and other devices;
- Essentially, an independent power circuit is required, and power sharing with other electronic devices is not recommended;
- Do not place the device in a location where it is difficult to remove the power plug or other removable plug.

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- The operator and manager of the device must fully understand the operating manual and it must be kept in close proximity to the equipment.
- To prevent safety-related accidents and ensure proper maintenance, the safety signs, instructions, and routine inspection chart included with the shock absorber should be posted in a clearly visible location near the shock absorber.
- Since there is a fan inside the UVT machine to circulate air, you should remove curtains or any other objects that may block the air flow near the UVT machine.
- When operating the UVT device, do not use mobile phones, radios, mobile radio transmitters, wireless toys, etc. in close proximity to the UVT device.
- During operation of the UVT device, the patient is prohibited from drinking drinks, water, etc., which may affect the operation of the UVT device.
- Do not separate the shock wave stimulation handle from the applicator while the shock wave machine is operating, as this may result in injury due to air pressure release.
- If transmission gel* is not applied sufficiently or does not make contact with the skin, patients may experience pain because an air bubble is formed and a shock wave is generated. Therefore, it is necessary to apply a sufficient amount of gel and press the shock wave stimulation handle firmly onto the patient to prevent pain.

* As a replacement for CARE SONIC transmission gel, the following gels can be used:

- Ultrasound gel " Mediagel " with Aloe juice Vera according to TU 21.20.23-001-76063983-2019, manufactured by Geltek -Medika LLC, Russia, registration certificate RZN 2021/15126 dated 08/19/2021.

- Medical gel for physiotherapy " Akugel-Rapan " according to TU 32.50.50-014-66242897-2021, produced by MediCraft LLC , Russia, registration certificate RZN 2023/20993 dated 08/29/2023.

- Shockwave may be subject to overpressure when used around the lungs, ears, and gastrointestinal tract.
- Do not apply shock waves directly to the heart, head, or injured area.
- The SWT device should be used for patients who can accurately express and convey pain during the procedure.
- Excessive pressure on the treatment area with the SWT stimulation handle may cause skin damage such as blistering and redness.
- The SWT machine should not be moved during use because the SWT stimulating handle may fall or the SWT machine may be damaged due to impact or movement of the SWT machine.

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Principle of action

A shock wave is a high-pressure wave that occurs when energy is suddenly released in a small space for a short period of time, stimulating damaged tissue and resulting in a therapeutic effect. Shock waves stimulate the early expression of angiogenesis-related factors, eNOS (endothelial nitric oxide synthase), VEGF (vascular endothelial growth factor) and PCNA (proliferating cell antinuclear antigen). It promotes neovascularization and affects blood supply and cell proliferation. As a result, it activates the surrounding tissues and bones to reduce pain and enhance natural healing effects. It also helps rebuild bones.

Shockwave is also known to have a therapeutic effect on chronic inflammation. Chronic inflammation is caused by the inflammatory response not being completely interrupted, which can cause damage to healthy tissue, resulting in chronic inflammation. Mast cells are the most important factor in the inflammatory process. The shock wave activates these cells, stopping the inflammatory response and aiding normal healing and repair processes. Shockwave also affects collagen, which is the most important factor in the process of repairing the upper musculoskeletal system and ligament. It also promotes the formation of procollagen, making the newly formed tendon firmer and firmer. In addition, shockwave is known to help treat pain by reducing the concentration of substance P, which is associated with pain signals.

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1 PURPOSE AND CONTRAINDICATIONS

Shock wave therapy (SWT) device, model SALUS - RSWT **Designed for** shock wave therapy. For more details, see "indications for use"

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Indications for use:

Neurology

Pain syndromes: dorsalgia (back), cervicgia (neck), trigeminalgia (face), headaches, cervical tension, greater occipital nerve compression syndrome (neck), piriformis syndrome (buttock), intercostal neuralgia (back), tunnel syndromes (upper and lower limbs, back, neck) .

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1. Do not use the device in conjunction with other electronic medical devices;
2. The following patients must be approved by a physician for the procedure:
patients with high body temperature (more than 37 °C), pregnant women, elderly people (over 85 years old) and children (under 3 years old);
a person with implants (eg, artificial hip joint, artificial pacemaker, implantable infusion port system, hearing aid).
3. It is prohibited to treat the following patients with the product:
blood disorders, clotting problems, or use of anticoagulants;
blood thinners (warfarin);
thrombosis;
oncological diseases;
pregnant women;
acute inflammation in the target area;
growth of cartilage tissue in children;
contact infectious patients (to avoid transmission of infection);
corticoid therapy for 6 weeks.

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
The shock wave therapy device consists of a housing , a power cord, a stimulating handle for shock wave therapy, and applicators for shock wave therapy.

The housing of the UVT device is equipped with a switch with a key, an emergency shutdown button, a touchscreen LCD display (LCD -10.2") with a joystick , a compartment for storing applicators, rollers for movement with and without brakes, connectors for connecting the stimulating handle for UVT and the power cord, as well as pre-installed software (software).

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The stimulating handle for UVT consists of a housing, a hammer, a guide and a rubber insert, and is also equipped with a cable for connecting to the body of the UVT device.

Applicators for shock wave therapy are manufactured in the following sizes: 6mm, 10mm, 15mm, 20mm and 35mm.

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
Figure 1. Front of the device

Stimulating handle for
shockwave therapy

Stimulating handle for
shockwave therapy


Control panel: 7-segment display, control buttons and
their movement control

Shock wave induction unit connector in channel 2

Shock wave induction unit connector in channel 1

Figure 2. Rear side of the device

AC receiver

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Power Connector

Power Switch

Figure 3 . Stimulating handle for UVT


Figure 4 . Control Panel Launch Screen


Figure 5. Power cord connector with bracket holder


Figure 6. Applicator appearance and size

6mm	10mm	15mm	20mm	35mm(option)
				

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Note: Firmware and software perform different functions together in one product. Firmware software is designed to control the hardware - controls the output in accordance with the input

from the software), the software is designed to control only what is shown on the display: therapy time, buttons.

⚠arning

- **Software updates should only be performed by authorized personnel. The software is encrypted and cannot be changed independently.**
-

EQUIPMENT:

Shock wave therapy (SWT) device, model SALUS - RSWT - DUET, consisting of:

- 1) UVT device, model SALUS-RSWT-DUET – 1 piece;
- 2) Stimulating handle for UVT – 2 pcs;
- 3) Power cord – 1 piece;
- 4) Applicator 6 mm – 1 piece;
- 5) Applicator 10 mm – 1 piece;
- 6) Applicator 15 mm – 1 piece;
- 7) Applicator 20 mm – 1 piece;
- 8) Applicator 35 mm – 1 piece;
- 9) Transmission gel CARE SONIC , 250 ml – 1 bottle;
- 10) Instructions For Use – 1 piece ;
- 11) Appendix to the Instructions For Use – 1 pc.

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Parameter name	Meaning, norm
1	2
Overall dimensions of the UVT device (L×W×H), mm • SALUS-RSWT-DUET	993 X 467 X 501
Roller diameter with brake, mm • SALUS-RSWT-DUET	75
Number of rollers with brake, pcs. • SALUS-RSWT-DUET	2
Roller diameter without brake, mm • SALUS-RSWT-DUET	75
Number of rollers without brake, pcs. • SALUS-RSWT-DUET	2
Force required to move the product, N • SALUS-RSWT-DUET	200
Weight of UVT device, kg • SALUS-RSWT-DUET	63
Dimensions of the stimulating handle for UVT (H×L×W), mm • SALUS-RSWT-DUET	33.8×195.5×86.5
Weight of the stimulating handle for UVT, kg • SALUS-RSWT-DUET	0.6
Cable length of the stimulating handle for UVT, mm • SALUS-RSWT-DUET	2000
Dimensions of the stimulating handle holder for UVT (H×D×W), mm • SALUS-RSWT-DUET	49 × 77 × 106.4
Number of stimulating handles for UVT, pcs. • SALUS-RSWT-DUET	2
Power cord length, mm • SALUS-RSWT-DUET	1500
Power cord weight, g • SALUS-RSWT-DUET	130
Number of operating modes, pcs • SALUS-RSWT-DUET	1 (shock wave therapy)
Time of readiness for work after switching on, min., no more • SALUS-RSWT-DUET	10
Maximum continuous operation time, min., not less Working mode	

SALUS-RSWT-DUET	Long-term with UVT
Rated supply voltage, V	
SALUS-RSWT-DUET	230
Power consumption, VA, no more	
SALUS-RSWT-DUET	1600
Rated value of supply frequency, Hz	
SALUS-RSWT-DUET	50
Display type	
SALUS-RSWT-DUET	built-in
Joystick dimensions (Diameter×Height), mm	
SALUS-RSWT-DUET	34×14
Dimensions of the emergency stop button (Diameter×Height), mm	
SALUS-RSWT-DUET	40×28
Dimensions of the base of the switch for the key (Diameter × Height), mm	
SALUS-RSWT-DUET	30×11
Key switch dimensions (Diameter×Height), mm	
SALUS-RSWT-DUET	21×14
Display	
SALUS-RSWT-DUET	10.2" LCD, touch screen
Display resolution	
SALUS-RSWT-DUET	800x480 pixels
Overall dimensions of fuses (Ø×D), mm	
SALUS-RSWT-DUET	6×30
Fuse current, A	
SALUS-RSWT-DUET	20
Fuse voltage, V	
SALUS-RSWT-DUET	250

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Table 3 Characteristics of applicators

Applicator size	Dimensions (height×diameter), mm	Impact surface diameter, mm	Weight, g
1	2		3
Applicator 6 mm	40.0×31.0	6	31
Applicator 10 mm	47.1×31.0	10	42.5
Applicator 15 mm	38.7×31.0	15	57
Applicator 20 mm	37.1×36.2	20	68
Applicator 35 mm	38.1×44.0	35	100.5

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The UVT stimulation handle is designed to be cleanable using the included cleaning tools.

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The body and elements of the SALUS-RSWT-DUET body, the stimulating handle for UVT are Cold stainless steel Rolled Coil , manufactured by HYUNDAI STEEL, Korea, ABS plastic, grade AF365, manufactured by LG Chemical , Korea.

Power cord, stimulating handle cable for UVT - Polyvinyl chloride brand CBL-REWT-GUN, manufactured by K TECHNOLOGY, Korea;

Applicators - 1) stainless steel SUS630 grade 630 H1100, manufactured by SeAH , Korea, 2) black polyacetal grade POM-C, manufactured by Namil Enpla CO., LTD., Korea;

Transmission gel CARE SONIC, 250 ml - Composition: ethanol, triethanolamine , glycerin, polyacrylic acid, water. Color: transparent, country of origin: Republic of Korea.

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In the manufacture of UHT devices, epoxy- polyester dyes of the EW096K brand, produced by AkzoNobel , Korea are used: white, black, gray, yellow, blue.

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The operating mode parameters correspond to Table 4.


Table 4

Parameter name	Meaning, norm
1	2
<i>Shock wave therapy parameters:</i>	
Operating modes SALUS-RSWT-DUET	Manual with the ability to store 9 programs in the device's memory with adjusted parameters of shock wave therapy. Adjustable parameters: number of beats, frequency, power.
Frequency range, Hz SALUS-RSWT-DUET with discreteness, Hz and permissible deviation, %	1-21 1 ± 20
Pressure range on the stimulating handle for UVT SALUS-RSWT-DUET , bar with discreteness, bar and permissible deviation, %	1-6 1 ± 20
Range of number of strokes, strokes SALUS-RSWT-DUET with discreteness, beats	100-5000 100
Standard sizes of used applicators, mm SALUS-RSWT-DUET	6, 10, 15, 20, 35

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Disease/location of affected areas	Processing frequency (Hz)	Pressure (intensity), bar	Mode
Plantar fasciitis (Plantar fasciitis)	13	2.5	1
Fascia pain (Musculofascial pain)	17	2.5	2

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
Table 6. Operating modes of the shock wave therapy (SWT) device, model: SALUS-RSWT-DUET

Disease/location of affected areas	Processing frequency (Hz)	Number of hits	Pressure (intensity), bar	Mode
Dorsalgia , cervicalgia , greater occipital nerve compression syndrome, piriformis syndrome, flat feet, scoliosis, inactivity atrophy syndrome, increased muscle tone, dystonia, pain, fibromyogia , brachial syndrome, treatment of fibromyalgia ,	5	2000-4000	1-3	Manual
Trigeminalgia , intercostal neuralgia, tunnel syndromes, polyneuropathy , intervertebral hernia,	5	2000-3000	1-2	Manual

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23

hyposthesia , rehabilitation after neurosurgical operations, rehabilitation after maxillofacial surgery, thick fiber pain, fine-fiber pain neuropathic pain, phantom pain, treatment of chronic fatigue				
Facial nerve paresis peroneal nerve palsy, radial nerve palsy, median nerve palsy, ulnar nerve paresis, recurrent nerve paresis, plexopathy	5	2000-4000	1-3	Manual
spinal injury, consequences of a stroke of the brain and spinal cord, consequences of operations on the brain and spinal cord	5	2000-4000	1-2	Manual
Arthritis, arthrosis, rehabilitation after arthroscopic operations, rehabilitation after injuries, epicondylitis , sacroiliac joint syndrome, pathology of the temporomandibular joint, ligamentites , traction degeneration of the plantar fascial band at the beginning of the medial tubercle of the calcaneus, glenohumeral periarthritis	5-10	2000-4000	1-6	Manual
Lymphedema, decongestant therapy	10	2000-3000	1-2	Manual
Overactive bladder, erectile dysfunction, premature ejaculation, acute and chronic pelvic pain syndrome,	5	2000	1-3	Manual

rehabilitation after prostate surgery, hypotension of the bladder, detrusor -sphincter dyssynergia , hyposthesia of the glans penis, pathology of the pudendal nerve, stress urinary incontinence, rehabilitation of weak pelvic muscles, restoration of neuromuscular control in the treatment of urinary incontinence.				
Irritable bowel syndrome, dysphagia	5Hz	2000	1-2	Manual
Physiotherapy for fractures, treatment of osteoporosis	5	2000-4000	1-2	Manual

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Periodically clean the outside of the shock wave machine, shock handle and touch screen with a soft cloth moistened with 70% isopropyl alcohol. Do not use abrasives, varnishes, thinners, ethylene or oxides as these materials may cause permanent damage to the UVT device;

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10.1 Safety requirements

Electrical safety precautions for the SALUS-RSWT-DUET model:

Check the adequacy of the power supply (230 VAC).

Check that all parts (power supply or peripherals) are connected correctly to the device.

Check ground continuity.

The device should be installed away from generators, X-ray machines, radio broadcast equipment, passing electrical wires, etc. to avoid possible noise during use.

Should use independent power supply and stable maintenance, it is not recommended to use the power supply with other electronic devices.

	Attention	The device should be installed away from generators, X-ray machines, radio broadcast equipment, passing electrical wires, etc. to avoid possible noise during use. Should use independent power supply and stable maintenance, it is not recommended to use the power supply with other electronic devices. When installing or operating the device, a specialist or authorized person must carry out periodic safety checks. For electrical safety reasons, the device must be connected to a grounded power source. Do not use the device in conjunction with other electronic medical devices. Do not use the device in areas where there are anesthetics or solvents. Do not install the device near substances such as oil, medicine, etc., which may pose a hazard to its operation. Do not use mobile phones, radios, mobile radio transmitters, wireless toys, etc. while operating the device. The patient should not consume drinks, water, etc., which may affect the operation of the device. If, while using the device, abnormalities in the body other than expected are observed, you should immediately stop using the device and consult a doctor.
---	------------------	--

Only special personnel duly authorized by the Manufacturer can repair and install the device; unauthorized disassembly or assembly of the device by the user is strictly prohibited.

If the user connects to the device other devices not specified in the Instructions for Use, this should be reported to the manufacturer or distributor authorized to distribute the device

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The outer sheath of the power cord of the UVT device, the cable of the stimulating handle for UVT should not peel off, the internal wires should not be exposed and should not be damaged by external shocks;

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12.4 Replacing consumable components of the stimulating handle for UVT

the hammer, guide and rubber insert of the stimulating handle for shock wave after every 1,500,000 - 3,000,000 (3 bar, 15 Hz) discharges of use, but the replacement period may vary depending on operating conditions.

12.5 To ensure and guarantee safe use, be sure to have the UHT device, including internal components and output voltage, inspected by personnel authorized by the manufacturer once a year.

12.6 Repairs to the UVT device may only be carried out by personnel authorized by the manufacturer.

12.7 In case of post-warranty and out-of-warranty repair of the UVT device, contact the authorized representative of the manufacturer in the Russian Federation: Limited Liability Company " Infomed-Neuro " (LLC " Infomed-Neuro ") at the address 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room I, fl . 1, room 50; phone +7 495 645 47 00; email davletov @ im - neuro . ru .

13 DISPOSAL INFORMATION

Disposal of the UVT device must be carried out in accordance with SanPin 2.1.3684-21 for class A waste, including applicators after preliminary disinfection, belongs to class A - epidemiological safe waste.

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
Authorized representative of the manufacturer in the Russian Federation: Limited Liability Company " Infomed-Neuro " (LLC " Infomed-Neuro ") at the address 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl . 1, room 50; phone +7 495 645 47 00; email davletov @ im - neuro . ru .

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Organization seller _____

Seller's signature _____

The product has been checked and has no damage.

I have read and agree with the operating rules.

Buyer's signature _____

The warranty period is 12 months from the date of sale from the manufacturer's warehouse.

For warranty questions, please contact: Limited Liability Company " Infomed-Neuro " (LLC " Infomed-Neuro ") at the address: Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl . 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .

In case of post-warranty and out-of-warranty repair of the product, contact: Limited Liability Company " Infomed-Neuro " (LLC " Infomed-Neuro ") at the address: 1 Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl . 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product/article must exactly match the number on the coupon. All columns of the coupon (except for the columns for the guarantee center) must be completed.

If the coupon is filled out inaccurately or incompletely, the product warranty is considered lost.

The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.

The product should not bear traces of mechanical or thermal influences, as well as traces of exposure to various chemical agents. Operation of the product is permitted only in accordance with the instruction manual for the product.

WARRANTY NOTE

DATE	DESCRIPTION	WARRANTY CENTER

15 INFORMATION ABOUT COMPLAINTS

In case of warranty and out-of-warranty repair of the product, please contact: Limited Liability Company " Infomed-Neuro " (LLC " Infomed-Neuro ") at the address: Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl . 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .

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SALUS-RSWT-DUET is intended for use in the electromagnetic environment described below. Purchasers or users of SALUS-RSWT-DUET must ensure that it is used under such conditions.

Emissions test	Correspondence	Electromagnetic Environment - Guidance
Radio interference GOST CISPR 11-2017	Group 1	SALUS-RSWT-DUET uses radio frequency energy only For their internal functions . That's why his radio emission Very low And hardly whether will cause any interference V located near electronic equipment
Radio interference GOST CISPR 11-2017	Class A	SALUS-RSWT-DUET suitable For use in everyone institutions except residential , and Maybe be used V household institutions And V those who directly connected To public low voltage networks power supply buildings used For household goals and self-care.
Harmonic radiation GOST IEC 61000-3-2-2021	Not applicable	
Voltage fluctuations/ flickering radiation GOST IEC 61000-3-3-2015	Not applicable	

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SALUS-RSWT-DUET is intended for use in the electromagnetic environment described below. Purchasers or users of SALUS-RSWT-DUET must ensure that it is used under such conditions.

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SALUS-RSWT-DUET is intended for use in the electromagnetic environment described below. Purchasers or users of SALUS-RSWT-DUET must ensure that it is used under such conditions.

Noise immunity test	Test Level IEC 60601	Compliance Level	Electromagnetic Environment - Guidance
Conductive radio frequency electromagnetic field GOST IEC 61000-4-6-2022	3 V r.m.s. 150 kHz – 80 MHz	3 V r.m.s.	Portable or mobile communications equipment should not be used on any cable section closer than this separation distance calculated

Resistance to RF fields GOST IEC 61000-4-3- 2016	3 V/m 80 MHz – 2.5 GHz	3 V/m	from the equation applicable to the transmitter frequency. 80 MHz – 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz – 2.5 GHz Where P is the maximum rated output power of the transmitter in watts (W) as specified by the transmitter manufacturer, and d is the recommended distance in meters (m). The field strength of a fixed RF transmitter, as determined from an electromagnetic field survey, must be less than the compliance level for each frequency range. B Interference may occur in the vicinity of equipment marked with the following symbol: 
--	---------------------------	-------	---

NOTE 1: At 80 MHz and 800 MHz applies more high frequency range .

NOTE 2: These recommendations can apply Not in everyone situations On spreading electromagnetic radiation influences absorption And reflection from structures , objects And people

^a Field strengths from fixed transmitters such as radio (cellular/cordless) telephones and base stations for land mobile radios, amateur radio, AM and FM radio and television broadcasts cannot be theoretically predicted with accuracy. To assess the electromagnetic environment created by a stationary radio transmitter, an electromagnetic field survey should be considered. If the measured field strength at the location where SALUS-RSWT-DUET is used exceeds applicable RF compliance level, monitoring is required to ensure that SALUS-RSWT-DUET is operating properly. If abnormalities are observed, additional measures such as orienting or adjusting the position of the SALUS-RSWT-DUET may be required.

b. The field strength must be less than [V1] V/m in the frequency range from 150 kHz to 80 MHz.

Recommended distance between portable and mobile radio communication devices and SALUS-RSWT-DUET

SALUS-RSWT-DUET is intended for use in an electromagnetic environment in which radiated radio frequency interference is controlled. The purchaser or user of SALUS-RSWT-DUET can help prevent electromagnetic interference by maintaining a minimum separation distance between SALUS-RSWT-DUET portable and mobile RF communications equipment (transmitters) as recommended below, consistent with the maximum communications output power.			
Rated maximum output power of the transmitter, W	Separation distance depending on transmitter frequency		
	150 kHz - 80 MHz	80 MHz - 800 MHz	800MHz - 2.5GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ Where $V_1 = 3$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ Where $E_1 = 3$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ Where $E_1 = 3$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

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NOTE 2: These guidelines may not apply in all situations. The propagation of electromagnetic radiation is affected by absorption and reflection from structures, objects and people. The manufacturer assumes no responsibility for interference caused by the use of a non-recommended connecting cable or by unauthorized modification or alteration of the device. Unauthorized modifications or alterations to the UVT equipment may void the user's authority to operate the UVT equipment.

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JANG MYEONG SU attorney-in-fact of
REMEDI CO., LTD.
President / Geun-Yong, Lee

appeared before me and admitted said principal's
subscription to the attached
DECLARATION.

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This is hereby attested on this
5th day of Jun. 2024 at this office.

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Daehan Law Firm & Notary Office Inc.
Belong to Seoul Central District
Prosecutors' Office
DaeGyung B/D 4F, 26-1 Sangdo-ro,
Dongjak-gu, Seoul, Republic of Korea

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Hyuk Jung

KWON HYUK JUNG

This office has been authorized by the
Minister of Justice, the Republic of
Korea, to act as Notary Public since
January 6, 2010 under Law No. 11823

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FAX : (02) 3477-0957

DECLARATION

We (applicant) : REMED CO., LTD.
do hereby solemnly and sincerely declare

1. That the attached document(s) Appendix to the Instructions For Use
SALUS-FSWT


is(are) true and correct.

is(are) true translation for the text(s) originally written in the Korean
language conscientiously believing the same to be true and correct.

2024. 06. 04.

Signature

APPROVE


CEO
«REMED Co., Ltd.»
Eun Hyun Koh, Sang Yong Lee

4th June 2024

Appendix to the Instructions For Use

Shock wave therapy (SWT) machine model SALUS-FSWT

production: REMED Co., Ltd., Republic of Korea

Appendix to the Operating Manual

Shock wave therapy (SWT) device, model SALUS-FSWT
(hereinafter referred to as UVT device(s), product(s), device(s), SALUS-FSWT)

Before using this device, please read the instruction manual and this manual appendix carefully and keep it for future reference!

SAFETY INSTRUCTIONS

ELECTRICAL SAFETY PRECAUTIONS:

- The rated power supply voltage for the UVT device is 220 V, alternating current;
- Make sure that all connecting parts (power cord or component cables) are correctly connected to the UVT device;
- Make sure that the UVT device is completely grounded;
- Repairs, extensions and installations of the UVT apparatus must not be carried out by anyone other than specialized personnel authorized by the manufacturer. Arbitrary disassembly/assembly of the UVT device by the user is strictly prohibited;
- Before connecting other equipment not specified in the instruction manual, be sure to consult the manufacturer or the manufacturer's authorized representative;
- To avoid electrical interference during use, the UVT device should be installed at a considerable distance from the power source, any generator, X-ray equipment, transmitting device, mobile electrical wire, etc.;
- Be careful because of possible electromagnetic or other interference between the UVT device and other devices;
- Essentially, an independent power circuit is required, and power sharing with other electronic devices is not recommended;
- Do not place the device in a location where it is difficult to remove the power plug or other removable plug.

BEFORE USING THE DEVICE, READ THE INSTRUCTION MANUAL AND THIS APPENDIX

WARNINGS AND PRECAUTIONS DURING OPERATION:

- The installation and re-installation of the UVT device must always be carried out by specialized personnel authorized by the manufacturer.
- Attention! Make sure that the patient does not have any metal objects during the procedure. Failure to comply with this condition during the procedure can lead to serious skin burns!
- The installed UVT device is subject to regular safety checks by specialized personnel authorized by the manufacturer.
- Repair or installation of the UHT device must only be carried out by specialized personnel authorized by the manufacturer, and therefore any arbitrary disassembly/assembly of the equipment by the user is strictly prohibited.
- For electrical safety reasons, the UVT device must always be connected to a reliably grounded power source for operation.

- The operator and manager of the device must fully understand the operating manual and it must be kept in close proximity to the equipment.
 - To prevent safety-related accidents and ensure proper maintenance, the safety signs, instructions, and routine inspection chart included with the shock absorber should be posted in a clearly visible location near the shock absorber.
 - Since there is a fan inside the UVT machine to circulate air, you should remove curtains or any other objects that may block the air flow near the UVT machine.
 - When operating the UVT device, do not use mobile phones, radios, mobile radio transmitters, wireless toys, etc. in close proximity to the UVT device.
 - During operation of the UVT device, the patient is prohibited from drinking drinks, water, etc., which may affect the operation of the UVT device.
 - Do not separate the piezo tip from the gel pad while the UVT device is in operation, as this may result in injury due to air pressure release.
 - If transmission gel* is not applied sufficiently or does not make contact with the skin, patients may experience pain because an air bubble is formed and a shock wave is generated. Therefore, it is necessary to apply a sufficient amount of gel and press the piezo tip firmly onto the patient to prevent pain.
- * As a replacement for CARE SONIC transmission gel, the following gels can be used:
- Ultrasound gel " Mediagel " with Aloe juice Vera according to TU 21.20.23-001-76063983-2019, manufactured by Geltek -Medika LLC, Russia, registration certificate RZN 2021/15126 dated 08/19/2021.
 - Medical gel for physiotherapy " Akugel-Rapan " according to TU 32.50.50-014-66242897-2021, produced by MediCraft LLC , Russia, registration certificate RZN 2023/20993 dated 08/29/2023.
- Shockwave may be subject to overpressure when used around the lungs, ears, and gastrointestinal tract.
 - Do not apply shock waves directly to the heart, head, or injured area.
 - The SWT device should be used for patients who can accurately express and convey pain during the procedure.
 - Excessive pressure on the treatment area with the piezo handpiece may cause skin damage such as blisters and redness.
 - Do not move the shock treatment device during use because the piezo tip may fall or the shock treatment device may be damaged due to impact or movement of the shock treatment device.

Do not use the UVT device in any of the following locations:

- in places where the ambient temperature or humidity does not meet the requirements;
- in areas where corrosive or flammable gases may be present;
- in places where there may be strong shock or vibration;
- in places directly exposed to wind and rain;
- in places where dust accumulates.

When operating the UVT device in the following locations, provide appropriate and sufficient protective measures:

- in areas exposed to static electricity or any other interference;
- in places exposed to intense electric or magnetic fields;
- in places located near power lines.

Operating principle

A shock wave is a high-pressure wave that occurs when energy is suddenly released in a small space for a short period of time, stimulating damaged tissue and resulting in a therapeutic effect. shock waves stimulate the early expression of angiogenesis- related factors , eNOS (endothelial nitric oxide synthase), VEGF (vascular endothelial growth factor) and PCNA (proliferating cell antinuclear antigen). It promotes neovascularization and affects blood supply and cell proliferation. As a result, it activates the surrounding tissues and bones to reduce pain and enhance natural healing effects. It also helps rebuild bones.

Shockwave is also known to have a therapeutic effect on chronic inflammation. Chronic inflammation is caused by the inflammatory response not being completely interrupted, which can cause damage to healthy tissue, resulting in chronic inflammation. Mast cells are the most important factor in the inflammatory process. The shock wave activates these cells, stopping the inflammatory response and aiding normal healing and repair processes. Shockwave also affects collagen, which is the most important factor in the process of repairing the upper musculoskeletal system and ligament. It also promotes the formation of procollagen , making the newly formed tendon firmer and firmer. In addition, shockwave is known to help treat pain by reducing the concentration of substance P, which is associated with pain signals.

INTRODUCTION

The main purpose of the manual and this instruction manual supplement is to provide the user with instructions for operating the UVT device, as well as information regarding its repair and maintenance. The manual provides the performance characteristics of the UVT device, its installation, operation and maintenance, as well as safety information. Before use, the user must carefully read the manual and this appendix to the operating manual in order to operate the UVT device correctly and accurately, achieving the specified safety and performance standards.

1 PURPOSE AND CONTRAINDICATIONS

Shock wave therapy (SWT) device, model SALUS - FSWT **Designed for** shock wave therapy. For more details, see "indications for use"

Conditions of use: UVT devices are suitable for professional use in hospitals and clinics by physiotherapists and doctors who have documents authorized to use this method.

Area of application : neurology, traumatology, urology and gynecology, gastroenterology, cosmetology, maxillofacial surgery.

Indications for use:

Neurology

Pain syndromes: dorsalgia (back), cervicgia (neck), trigeminalgia (face), headaches, cervical tension, greater occipital nerve compression syndrome (neck), piriformis syndrome (buttock), intercostal neuralgia (back), tunnel syndromes (upper and lower limbs, back , neck) .

Peripheral paresis: facial nerve palsy (face), peroneal nerve palsy (lower limb), radial nerve palsy (upper limb), median nerve palsy (upper limb), ulnar nerve palsy (upper limb), recurrent nerve palsy (neck).

Central paresis: spinal injury (any part of the body except the head (scalp)), consequences of a stroke of the brain and spinal cord (any part of the body except the head (scalp)), consequences of operations on the brain and spinal cord (any part of the body except the head (scalp) heads)).

Plexopathies (any part of the body except the head (scalp))

Complex treatment of the consequences of intervertebral hernia (any part of the body except the head (scalp))

Hyposthesia (any part of the body except the head (scalp))

Rehabilitation after neurosurgical operations (any part of the body except the head (scalp))

Traumatology

Joint diseases:

Arthrosis (any part of the body except the head (scalp)), arthritis (any part of the body except the head (scalp)) .

Rehabilitation: after arthroscopic operations (any part of the body except the head (scalp)), after injuries (any part of the body except the head (scalp)).

Pain syndromes: impingement syndrome (shoulder, upper limb, neck), myofascial syndrome (any part of the body except the head (scalp)), epicondylitis (upper limb shoulder, neck), plantar fasciitis (lower limb, back), tibial tract syndrome (lower limb, back), sacroiliac joint syndrome (pelvis, back, lower limb), trochanteritis (pelvis, back, lower limb), pathology of the temporomandibular joint (face) , pathology of the Achilles tendon (lower limb, back), ligamentitis (any part of the body except the head (scalp)) .

Movement disorders:

Flat feet (lower limb), scoliosis (back), dystonia (any part of the body except the head (scalp)), disuse atrophy syndrome (any part of the body except the head (scalp)).

Lymphatic edema .

Urology and gynecology

Rehabilitation after prostate surgery, erectile dysfunction , acute and chronic pelvic pain syndrome, **etrusor** - sphincter dyssynergia , stress urinary incontinence, Overactive bladder (pelvis, lower limb, back), Premature ejaculation (pelvis, lower limb, back), Bladder hypotension (pelvis, lower limb, back), Hypoesthesia of the glans penis (pelvis, lower limb, back), Pathology of the pudendal nerve (pelvis, lower limb, back)

Gastroenterology

- Irritable bowel syndrome (stomach, back)
- Dysphagia (neck, back)

Cosmetology

Increased muscle tone: abdomen, gluteal region, shoulder, face (any part of the body except the head (scalp)).

Decongestant therapy (any part of the body except the head (scalp)).

Maxillofacial surgery

Rehabilitation after maxillofacial surgery (face)

Other diseases

Providing a completely non-invasive high-intensity pulsed electromagnetic effect on the pelvic floor muscles with the aim of rehabilitating weak pelvic muscles and restoring neuromuscular control in the treatment of urinary incontinence in humans (pelvis, back);

Shoulder syndrome (upper limb);

Humeroscapular periarthritis (upper limb);

Physiotherapy for fractures (any part of the body except the head);

Chronic proximal plantar fasciitis in patients aged 18 years or older with symptoms for 6 months or more, plantar fasciitis (lower limb);

Traction degeneration of the plantar fascial band at the beginning of the medial tubercle of the calcaneus (lower limb);

Repeated muscle contractions by passing electrical currents through electrodes without contact with the affected area of the body;

Relief of pain (myofascial pain, thick-fiber pain , thin-fiber pain , neuropathic pain , phantom pain, fibromyogia) (any part of the body except the head);

Treatment of osteoporosis (any part of the body except the head);

Treatment of fibromyalgia , arthritis, chronic fatigue (any part of the body except the head).

Providing a completely non-invasive, high-intensity pulsed electromagnetic effect on the pelvic floor muscles with the aim of rehabilitating weak pelvic muscles and restoring neuromuscular control in the treatment of urinary incontinence in humans (pelvis).

Pain in the ankle and calf muscles (lower limb), patellar tendinitis (lower limb), hip pain (lower limb), patellofemoral pain syndrome (lower limb), achillodynia (lower limb), hypertrophic tendinitis (any part of the body other than the head (scalp)), muscle tendon calcification (any part body except the head (scalp)), pain in the supraspinatus muscle (upper limb, shoulder), pain in the forearm and wrist (upper limb).

Contraindications for use:

Contraindications for use in shock wave therapy are as follows:

1. Do not use the device in conjunction with other electronic medical devices;
2. The following patients must be approved by a physician for the procedure:
patients with high body temperature (more than 37 °C), pregnant women, elderly people (over 85 years old) and children (under 3 years old);
a person with implants (eg, artificial hip joint, artificial pacemaker, implantable infusion port system, hearing aid).
3. It is prohibited to treat the following patients with the product:
blood disorders, clotting problems, or use of anticoagulants;
blood thinning drugs (warfarin);
thrombosis;
oncological diseases;
pregnant women;
acute inflammation in the target area;
growth of cartilage tissue in children;
contact infectious patients (to avoid transmission of infection);
corticoid therapy for 6 weeks.

Attention! Shock waves should not be applied over air-filled tissue (lungs), nor should they be used near large nerves, blood vessels, vertebrae, or on the head.

Possible complications:

If any abnormalities are found, such as skin redness, pain, hematoma with petechiae, hematoma, bleeding, swelling, you should immediately stop using the device and consult a doctor.

Possible side effects:

Side effects of shock wave therapy:

Swelling, redness, hematoma, pain.

2 GENERAL DESCRIPTION OF THE UVT DEVICE

The UVT device consists of a housing , a power cord, a piezo tip, a foot switch, and gel pads for shock wave therapy.

The case of the UVT device is equipped with a switch with a key, an emergency shutdown button, a touch LCD display (LCD -10.2") with a joystick , a compartment for storing gel pads, connectors for connecting a piezo tip and a power cord, as well as pre-installed software (software)


The display of the UVT device is built into the body and is located on the top panel.

Gel pads for shock wave therapy are manufactured in the following sizes: 5mm, 10mm, 15mm, 20mm, 25mm, 30mm, 35mm and 40mm.

3 DESIGN FEATURES AND APPEARANCE OF THE UVT DEVICE


DESIGN FEATURES AND APPEARANCE OF THE UVT DEVICE

Figure 1. Front of the device

No.	Name	Function
1	Tip connector	Connection of tip to main body
2	LCD touch screen	Screen display using touch screen
3	Encoder	Adjusting and setting stimulation data and intensity

Figure 2. Rear side of the device

No.	Name	Function
4.	Foot switch connector	Connecting a foot switch to the main unit
5.	Equipotential ground connector	Protecting people from electric shock
6.	Speaker	Voice output
7.	Ventilation hole	Ventilation for heat removal
8.	AC power input connector	AC power connector
9.	Switch	Device ON/OFF switch

Figure 3. Sides of the device




Figure 4. Piezo tip


No.	Name	Function
80.	Connector	Tip connector to main body
81.	Pen	Handle tip
82.	Shock wave generator	Shock wave generation
83.	Cable	Cable for connecting the handpiece to the main unit
89.	Gel pad for piezo tip	The penetration depth of the shock wave can be determined by the height of the gel pad

Figure 5. Foot switch

No.	Name	Function
84	Foot switch to control shock wave generation	Shock wave generation control
85	Foot switch to control downward adjustment	Zooming out of the selected item on the display
86	Foot switch to control up adjustment	Enlarge the selected item on the display
87	Foot switch connector	Connector for connecting the foot switch to the main body
88	Foot switch cable	Cable, foot switch connections to main unit

Figure 6. Control Panel (LCD Display)

No.	Function
31	Window displaying information about the currently selected mode
32	Output intensity display window
33	Window displaying information about the frequency of shock wave generation
34	Window displaying information about the frequency of shock wave occurrence
35	Button to change elements
36	Button to change elements
37	Button to start the selected mode
38	Settings button
39	Treatment time display window (a window that is calculated and displayed based on a given frequency and number of occurrences)
40	Remaining time (window for calculating and displaying the time required to display all cases)
41	Display window with selected tip
42	Electric handpiece selection button
48	Assist tool (touch and position change can be changed by MEMO) combined with 31 (mode information)
69	Zoom button

The device has pre-installed proprietary software (hereinafter referred to as the software). After turning on the product, the software automatically starts. The software life cycle diagram is presented in Appendix A.1. Information about the software corresponds to Table 1.

Table 1: Software Details

Software name	Software identification name	Software version number (identification number)	Digital software identifier (executable code checksum)	Algorithm for calculating the digital software identifier	Software release date
SALUS - FSWT Main software version	-	1.00.00	-	-	11/23/2018
SALUS - FSWT Slave version (FOCUS firmware)	-	1.00.00	-	-	11/23/2018

Note: Firmware and software perform different functions together in one product. Firmware software is designed to control the hardware - controls the output in accordance with the input from the software), the software is designed to control only what is shown on the display: therapy time, buttons.


Warning

- **Software updates should only be performed by authorized personnel. The software is encrypted and cannot be changed independently.**

EQUIPMENT:

Shock wave therapy (SWT) device, model SALUS - FSWT, consisting of:

- 1) UVT device, model SALUS-FSWT – 1 piece;
- 2) Piezo tip – 1 piece;
- 3) Foot switch – 1 piece;
- 4) Power cord – 1 piece;
- 5) Trolley (if necessary)
- 5) Gel pad 5 mm – 1 piece;
- 6) Gel pad 10 mm – 1 piece;

REMED
Rehabilitation Medical Company

-
- 7) Gel pad 15 mm – 1 piece;
 - 8) Gel pad 20 mm – 1 piece;
 - 9) Gel pad 25 mm – 1 piece;
 - 10) Gel pad 30 mm – 1 piece;
 - 11) Gel pad 35 mm – 1 piece;
 - 12) Gel pad 40 mm – 1 piece;
 - 13) Gel pad holder – 1 pc;
 - 14) Transmission gel CARE SONIC , 250 ml – 1 bottle;
 - 14) Instructions For Use – 1 piece;
 - 15) Appendix to the Instructions For Use – 1 pc.

4 TECHNICAL SPECIFICATIONS

The operational and technical characteristics of the UVT device correspond to Tables 2-3.
Table 2

Parameter name	Meaning, norm
1	2
Overall dimensions of the UVT device (L×W×H), mm • SALUS-FSWT	424x405x207
Roller diameter with brake, mm • cart	75
Number of rollers with brake, pcs. • cart	2
Roller diameter without brake, mm • cart	75
Number of rollers without brake, pcs. • cart	2
Force required to move the product, N • trolley with SALUS device - FSWT	115
Weight of UVT device, kg • SALUS-FSWT	13
Piezo tip dimensions (L×W×H), mm • SALUS-FSWT	116 x 116 x 151
Weight of piezo tip, kg • SALUS-FSWT	1
Piezo tip cable length, mm • SALUS-FSWT	2500
Foot switch dimensions, (L×W×H), mm • SALUS-FSWT	450 x 143 x 48
Foot switch weight, kg • SALUS - FSWT	2
Number of piezo tips, pcs • SALUS-FSWT	1
Power cord length, mm • SALUS-FSWT	1500
Power cord weight, g • SALUS-FSWT	130
Trolley dimensions, (L×W×H), mm SALUS-FSWT	625×470×770
Trolley weight, kg SALUS - FSWT	25
Dimensions of gel pad holder, (L×W), mm	

• SALUS-FSWT	350×240
Number of operating modes, pcs • SALUS-FSWT	(shock wave therapy)
Time of readiness for work after switching on, min., no more • SALUS-FSWT	1
Maximum continuous operation time, min., not less Working mode • SALUS-FSWT	Long-term with UVT
Rated supply voltage, V • SALUS-FSWT	220
Power consumption, VA, no more • SALUS-FSWT	300
Rated value of supply frequency, Hz • SALUS-FSWT	50
Display type • SALUS-FSWT	built-in
Joystick dimensions (Diameter×Height), mm • SALUS-FSWT	34×14
Dimensions of emergency stop button (Diameter×Height), mm • SALUS-FSWT	40×28
Dimensions of the base of the switch for the key (Diameter × Height), mm • SALUS-FSWT	30×11
Key switch dimensions (Diameter×Height), mm • SALUS-FSWT	21×14
Display • SALUS-FSWT	10.2" LCD, touch screen
Display resolution • SALUS-FSWT	800x480 pixels
Overall dimensions of fuses (Ø×D), mm • SALUS-FSWT	5×20
Fuse current, A • SALUS-FSWT	2
Fuse voltage, V • SALUS-FSWT	250

Note – Tolerances for all characteristics are $\pm 2.5\%$ unless otherwise stated.

Table 3 Characteristics of gel pillows

Gel pillow size	Diameter, mm	Height, mm	Weight, g
-----------------	--------------	------------	-----------

1	2	3	4
Gel pad 5 mm	5	6 2.8	219
Gel pad 10 mm	10	57.8	214
Gel pad 15 mm	15	52.8	207
Gel pad 20 mm	20	47.8	196
Gel pad 25 mm	25	37.8	158
Gel pad 30 mm	30	32.8	130
Gel pad 35 mm	35	23.2	99
Gel pad 40 mm	40	23.2	99

Note – Tolerances for all characteristics are $\pm 2.5\%$ unless otherwise stated.

The UVT device must remain operational at temperatures from $+10\text{ }^{\circ}\text{C}$ to $+40\text{ }^{\circ}\text{C}$, relative air humidity from 30 to 85% and atmospheric pressure from 700 to 1060 hPa.

The UVT device, packaged in a transport container, must be resistant to climatic influences, at temperatures from $-10\text{ }^{\circ}\text{C}$ to $+50\text{ }^{\circ}\text{C}$, relative air humidity from 0 to 85% and atmospheric pressure from 700 to 1060 hPa.

The UVT device is resistant to impacts caused by rough handling according to IEC 60601-1 and mechanical stress, with the following parameters:

- vibration loads – frequency range 10-55 Hz with a displacement amplitude of 0.15 mm;

The operating mode must be long.

Protection class against electric shock according to IEC 60601-1: Class I, applied part Type BF;

The UVT device must meet the requirements of electromagnetic compatibility (class A, group 1).

Protection against penetration of water and dust of the device - IP X0.

The outer surfaces of the UHT device parts are resistant to repeated cleaning and disinfection using a chemical method using 70% isopropyl alcohol.

The design of the UVT apparatus and its parts is repairable.

Service life - 5 years The criterion for the limiting state of the UVT apparatus is the impossibility or impracticality of restoration.

The output signal shape is a sinusoidal signal.

5 MATERIALS USED IN MANUFACTURE

UVT devices are made of the following materials:

The body and elements of the SALUS-FSWT body, trolley, piezo tip are cold-rolled steel EG COIL, manufactured by POSCO, Korea, ABS plastic AF365, manufactured by LG Chemical, Korea.

Foot switch – polycarbonate brand GN1002FC, manufactured by LG Chemical, Korea;

Power cord, piezo tip cable, foot switch cord cable - Polyvinyl chloride brand CBL-REWT-GUN, manufactured by K TECHNOLOGY, Korea;

Gel pads – Polyethylene terephthalate 3 T, produced by SK CHEMICALS, Korea;

Transmission gel CARE SONIC, 250 ml - Composition: ethanol, triethanolamine, glycerin, polyacrylic acid, water. Color: transparent, country of origin: Republic of Korea.

UVT devices do not contain drugs or components of animal or human origin.

In the manufacture of UHT devices, epoxy-polyester dyes of the EW096K brand, produced by AkzoNobel, Korea are used: white, black, gray, yellow, blue.

All color options are provided by mixing basic dyes.

Type and duration of contact with the human body: short-term contact (no more than 24 hours) with the surface of intact skin.

6 OPERATING MODE PARAMETERS

The operating mode parameters correspond to Table 4.

Table 4

Parameter name	Meaning, norm
1	2
<i>Shock wave therapy parameters:</i>	
Operating modes SALUS-FSWT	Manual with the ability to store 20 programs in the device's memory with adjusted shock wave therapy parameters (from M 1 to M 20) Adjustable parameters: number of shocks, frequency, pressure
Frequency range, Hz SALUS-FSWT	1-12
with discreteness, Hz	1
and permissible deviation, %	± 20
Energy flux density range at the piezo tip SALUS-FSWT , mJ/mm²	0-1
with discreteness, mJ/mm ²	0.05
and permissible deviation, %	± 20
Range of number of strokes, strokes SALUS-FSWT	100-10000
with discreteness, beats	100
Standard sizes of used gel pillows, mm SALUS-FSWT	5, 10, 15, 20, 25, 30, 35, 40

The preset (standard) operating modes of the UVT device are reflected in Table 5.

Table 5. Operating modes of the shock wave therapy (SWT) device, model: FSWT

Disease/location of affected areas	Processing frequency (Hz)	Number of hits	Intensity, % (energy flux density, mJ)	Mode
Cervical pain V neck)	5	2000	40 (0.4)	1
Shoulder pain	7	2000	50 (0.5)	2
Calcific tendinitis (calcification of muscle	7	2000	45 (0.45)	3

tendons)				
Supraspinatus tendon (Supraspinatus muscle)	6	2000	40 (0.4)	4
Epicondylitis (Ulnar)	5	1500	25 (0.25)	5
Epicondylitis (Radial) Lateral epicondylitis)	5	1500	25 (0.25)	6
Forearm treatment	8	1500	30 (0.3)	7
Wrist	5	1500	20 (0.2)	8
Thumb pain (Pain in the thumb)	6	1500	20 (0.2)	9
Lumbar pain	8	2000	55 (0.55)	10
Thigh treatment	10	1500	60 (0.6)	11
Trochanter pain (Trochanteritis)	8	1500	60 (0.6)	12
Hip pain (Pain in the hip joint)	12	2000	65 (0.65)	13
Patellar Tendinitis tendonitis)	4	2000	25 (0.25)	14
Calf muscles	6	1500	55 (0.55)	15
Tibia pain (Tibial pain)	8	1500	35 (0.35)	16
Ankle (Ankle)	6	1500	20 (0.2)	17
Achilles tendinitis (Achilles tendon)	5	2000	30 (0.3)	18
Plantar fasciitis (Plantar fasciitis)	5	1500	40 (0.4)	19
Myofascial pain (Myofascial pain)	8	2000	40 (0.4)	20

List of customizable (manual) operating modes of UVT devices:

Table 6. Operating modes of the shock wave therapy (SWT) device, model: FSWT

Disease/location of affected areas	Processing frequency (Hz)	Number of hits	Intensity , % (energy flux density, mJ)	Mode
Dorsalgia , cervicalgia , greater occipital nerve	5	2000-4000	5-50 (0.05-0.5)	Manual

compression syndrome, piriformis syndrome, flat feet, scoliosis, inactivity atrophy syndrome, increased muscle tone, dystonia, pain, fibromyogia , brachial syndrome, treatment of fibromyalgia ,				
Trigeminalgia , intercostal neuralgia, tunnel syndromes, polyneuropathy , intervertebral hernia, hyposthesia , rehabilitation after neurosurgical operations, rehabilitation after maxillofacial surgery, thick fiber pain, fine-fiber pain neuropathic pain, phantom pain, treatment of chronic fatigue	5	2000-3000	10-20 (0.1-0.2)	Manual
Facial nerve paresis peroneal nerve palsy, radial nerve palsy, median nerve palsy, ulnar nerve paresis, recurrent nerve paresis, plexopathy	5	2000-4000	5-50 (0.05-0.5)	Manual
spinal injury, consequences of a stroke of the brain and spinal cord, consequences of operations on the brain and spinal cord	5	2000-4000	10-20 (0.1-0.2)	Manual
Arthritis, arthrosis, rehabilitation after arthroscopic operations, rehabilitation after injuries, epicondylitis , sacroiliac joint syndrome,	5-10	2000-4000	5-100 (0.05-1)	Manual

pathology of the temporomandibular joint, ligamentites, traction degeneration of the plantar fascial band at the beginning of the medial tubercle of the calcaneus, glenohumeral peri arthritis				
Lymphedema, decongestant therapy	10	2000-3000	10-30 (0.1-0.3)	Manual
Overactive bladder, erectile dysfunction, premature ejaculation, acute and chronic pelvic pain syndrome, rehabilitation after prostate surgery, hypotension of the bladder, detrusor -sphincter dyssynergia, hyposthesia of the glans penis, pathology of the pudendal nerve, stress urinary incontinence, rehabilitation of weak pelvic muscles, restoration of neuromuscular control in the treatment of urinary incontinence.	5	2000	5-50 (0.05-0.5)	Manual
Irritable bowel syndrome, dysphagia	5	2000	5-30 (0.05-0.3)	Manual
Physiotherapy for fractures, treatment of osteoporosis	5	2000-4000	5-20 (0.05-0.2)	Manual

7 RULES FOR CARE, CLEANING AND DISINFECTION

The UVT device can be cleaned in various ways. Avoid any damage or contamination of the UVT device using the methods recommended below. Before cleaning the UVT device, be sure to turn it off. Cleaning the inside of the UVT device body is permitted only to personnel authorized by the manufacturer.

Periodically clean the external surface of the UVT device, piezo tip and touch screen with a soft cloth moistened with 70% isopropyl alcohol. Do not use abrasives, varnishes, thinners, ethylene or oxides as these materials may cause permanent damage to the UVT device;

Do not immerse any part of the UVT device in liquid or detergent. In addition, liquid must not be allowed to enter the UVT device or its components;

After each shock wave therapy treatment, clean the applicator with 70% isopropyl alcohol and a smooth cloth;

If you choose to use a disinfectant or method not mentioned above, it may damage the UHT machine. Please note that if the UVT device is damaged by unauthorized materials, no warranty will be provided.

It is necessary to regularly clean the applicators once a week, as the gel used during treatments may stick to the applicator after prolonged use and may affect the generation of the shock wave.

Use running water to clean applicators. If the gel is difficult to wash off, it is effective to immerse the applicator in water for 2-3 hours and then rinse. After washing, remove the water with a soft dry cloth and dry at room temperature for more than 1 hour.

To ensure safe use, the internal cleaning of the UHT device housing must be performed once a year by personnel authorized by the manufacturer.

8 MARKING


8.1 Marking applied to the label glued to the body of the UVT device:

- UVT device model;
- manufacturer's logo;
- serial number;
- date of manufacture;
- rated voltage;
- AC;
- rated frequency;
- maximum power consumption;
- working part type BF
- general warning sign;
- refer to the service instructions;
- do not dispose of with household waste;
- name and address of the manufacturer;
- "Do not push" symbol;
- "Do not lean" symbol.

8.2 Markings applied to the packaging of the device:

Symbol	Meaning	Description
	Do not stack.	Indicates that items cannot be stacked vertically, either due to the nature of the shipping packaging or the nature of the items themselves.
	Keep away from rain	Indicates that shipping packaging should be stored away from rain and in dry conditions.
	Fragile	To indicate that the contents of the shipping package are fragile and should be handled with care.
	This side up	To indicate the correct vertical position of the transport package.
	Recycling	To indicate the location of a trash can or container

8.3 Explanation of all marking symbols:

Symbol	Meaning	Description
	Attention	The phrase "Warning" is used to inform the user of the fact that there is a reason that may result in a fatal error if ignored.
	Carefully	The phrase "Caution" is used to inform the user that failure to follow safety instructions may result in injury, although there is no danger to life.
	Type BF	Working part type BF . Medical product with an increased degree of protection and an isolated working part
	Maintenance instructions	Where there is a "Maintenance Instructions" sign, please follow the requirements specified in the Maintenance Instructions.
	Serial number	Designation of the serial number of the medical device.
	Manufacturer	Designation of the manufacturer of the medical device.
	Production date	Designation of the date of manufacture of the medical device.
	Off (power)	Designation of the position of the power switch corresponding to disconnection from the network
	On (power)	Designation of the position of the power switch corresponding to connection to the network
	AC	Indication on the product plate about the need to connect the product through the appropriate connector only to AC sources
	Do not dispose of with household waste	Instructions on the need for disposal at separate collection points for waste recovery and recycling
	Don't push	Do not push the product while the parking brake is on or in the on state.
	Don't lean your elbows	You cannot lean your elbows on the product. This may result in equipment damage or personal injury.

9 PACKAGING

The device, packed in a plastic bag, is placed in shock-absorbing foam inserts. On top of the shock-absorbing insert, in a special recess, is placed a cardboard box containing a power cord, bracket, instruction manual in Russian, and keys. The top is closed with a plastic bag. Then the outer cardboard box is put on.

10 SAFETY AND ENVIRONMENTAL REQUIREMENTS

10.1 Safety requirements

Electrical safety precautions for model SALUS-FSWT:

Make sure the device's power line meets the requirements.

Make sure all connections (power lines or peripherals) are properly connected to the product.

Make sure the device is completely grounded. For proper connection, see the connection section.

To avoid electrical interference during use, install it away from generators, X-ray equipment, broadcast equipment, mobile devices and so on.

Separate power circuits and consistent maintenance are important in this chapter, and it is recommended that you do not turn off the power as you would with any other electronic device.

Repair, modification, or installation of the equipment is prohibited except by authorized personnel authorized by the manufacturer. Never disassemble or assemble the equipment yourself.

When connecting other equipment not specified in the Operating Instructions, you must notify the distributor who has the right to sell the company or its products.

In terms of safety, the UVT device must comply with the requirements of IEC 60601-1.

In terms of electromagnetic compatibility, the UVT device must comply with the requirements of IEC 60601-1-2 (Class A, Group 1).

Serviceability in accordance with IEC 60601-1-6 and IEC 62366.

Software class according to IEC 62304 – B.

10.2 Environmental protection

The materials used in the manufacture of the UVT device should not pose a danger to life, human health and the environment both during operation and after the end of its service life.

The main type of possible hazardous impact on the environment is pollution of atmospheric air in populated areas, soils and waters as a result of:

- emergency leaks (spills) of production materials;
- unorganized disposal of waste on the territory of the manufacturer or outside it;
 - - arbitrary dumping of waste in places not intended for this purpose.

11 TRANSPORTATION AND STORAGE

11.1 Transportation conditions

Transportation of the UVT device can be carried out by all types of covered vehicles and the rules for the transportation of goods that apply to each type of transport.

Transportation of the UVT apparatus by sea must be carried out in accordance with the "Safety Rules for the Maritime Transportation of General Cargoes". Type of shipment - containers and small shipments.

The UVT device, during transportation, is resistant to climatic influences at ambient temperatures from -10°C to $+50^{\circ}\text{C}$, relative air humidity from 0 to 85% and atmospheric pressure from 70 to 106 kPa.

The UVT device, during transportation, is resistant to mechanical stress, with the following parameters:

vibration loads – frequency range 10-55 Hz with a displacement amplitude of 0.15 mm.

During transportation it is prohibited:

Push the package with the UVT device;

Sit on the package with the UVT device;

Stack packages with UVT devices.

11.2 Storage conditions

The UVT device, during storage, is resistant to climatic influences at ambient temperatures from -10°C to $+50^{\circ}\text{C}$, relative air humidity of no more than 80% and atmospheric pressure from 70 to 106 kPa.

Storage conditions:

Keep away from water;

Keep away from direct sunlight;

Do not store near heaters;

Do not store in dusty places;

Avoid locations subject to strong shock or vibration, chemicals, or explosive gases.

12 MAINTENANCE AND REPAIR

12.1 Routine inspection of the UVT apparatus:

The outer sheath of the power cord of the UVT device, piezo tip should not peel off, the internal wires should not be exposed and should not be damaged by external shocks;

Clean the outside of the UVT device to ensure that there are no foreign objects;

The various parts attached to the UVT apparatus should not shake.

12.2 After long-term storage, be sure to check the UVT device before using it.

12.3 Removing condensation from the air compressor

Condensation should be removed after every 1,500,000 shocks of use. Otherwise, it may cause malfunction or unpleasant odor. Depending on the ambient humidity, more frequent condensate drainage may be required.

12.4 To ensure and guarantee safe use, be sure to have the UHT device, including internal components and output voltage, inspected by personnel authorized by the manufacturer once a year.

12.5 Repairs to the UVT device may only be carried out by personnel authorized by the manufacturer.

12.6 In the case of post-warranty and out-of-warranty repairs of the UVT device, contact the authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room I, fl. 1, room 50; phone +7 495 645 47 00; email davletov@im-neuro.ru.

13 DISPOSAL INFORMATION

Disposal of the UVT device must be carried out in accordance with SanPin 2.1.3684-21 for class A waste, including gel pads after preliminary disinfection, belongs to class A - epidemiological safe waste.

The UVT device is subject to disposal in the following cases:

end of service life;

confirmation of facts and circumstances that pose a threat to the life and health of health workers and indicate a failure to fulfill the intended purpose.

All packaging, including transport packaging, must also be recycled. Paper, wood, polyethylene and plastic must be disposed of separately in accordance with the rules for collection, accounting and disposal of packaging.

14 MANUFACTURER WARRANTY

The UVT device is guaranteed for a period of 12 months from the date of sale from the manufacturer's warehouse.

The guaranteed shelf life is 1 year if storage conditions are met.

The UVT device is guaranteed to be free from defects related to manufacturing and quality of production materials. Any defective part will be replaced as long as the device has not been tampered with and the device has been properly operated during this period. Make sure that any malfunction does not occur due to improper care of the device or non-compliance with the operating instructions. In the event of a malfunction, subject to compliance with operating requirements in accordance with the operating manual, the manufacturer carries out routine repairs or replacement of the UVT device.

The manufacturer or the manufacturer's authorized representative in the Russian Federation is not responsible for any personal/property damage caused by erroneous use, operation for purposes other than the intended purpose, and negligence in servicing the product.

Please note that if the UHT device is damaged by the use of unauthorized cleaning and disinfection materials, no warranty will be provided.

Authorized representative of the manufacturer in the Russian Federation: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl. 1, room 50; phone +7 495 645 47 00; email davletov @ im - neuro . ru .

WARRANTY CARD

Product category _____
 Manufacturer _____
 Execution option _____
 Serial number/art. _____
 Date of sale _____
 Organization seller _____
 Seller's signature _____
 The product has been checked and has no damage.
 I have read and agree with the operating rules.
 Buyer's signature _____

Stamp
about sale

The warranty period is 12 months from the date of sale from the manufacturer's warehouse.
 For warranty questions, please contact: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address: Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl. 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .
 In case of post-warranty and out-of-warranty repair of the product, contact: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address: 1 Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl. 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .

ATTENTION! Make sure that the warranty card is filled out correctly. The serial number of the product/article must exactly match the number on the coupon. All columns of the coupon (except for the columns for the guarantee center) must be completed.
If the coupon is filled out inaccurately or incompletely, the product warranty is considered lost.
The manufacturer's warranty is carried out in accordance with the Civil Code of the Russian Federation on a general basis.
The product should not bear traces of mechanical or thermal influences, as well as traces of exposure to various chemical agents. Operation of the product is permitted only in accordance with the instruction manual for the product.

WARRANTY NOTE

DATE	DESCRIPTION	WARRANTY CENTER

15 INFORMATION ABOUT COMPLAINTS

In case of warranty or out-of-warranty repair of the product, please contact: Limited Liability Company "Infomed-Neuro" (LLC "Infomed-Neuro") at the address: Russia, 125445, Moscow, Leningradskoye Shosse, building 69, building 1, room. I, fl. 1, room 50; phone +7 495 645 47 00; email kitayeva @ im - neuro . ru .

16 ELECTROMAGNETIC COMPATIBILITY

1. Electromagnetic Compatibility Specifications

Definition: EMC (electromagnetic compatibility) is the ability of equipment to avoid electromagnetic interference generated by other devices without causing similar electromagnetic interference to other equipment.

Warning

Because this device generates, uses, and radiates radio frequency energy, it may interfere with the operation of other medical and non-medical devices and wireless communications. It is characterized by the creation of electromagnetic interference in the air or cable environment. This product is designed to fully comply with EMC standards and provides reasonable protection against such interference.

Warning

UHT equipment requires special measures to ensure **ELECTROMAGNETIC COMPATIBILITY** and must be installed and put into service in accordance with the EMC-related information given below. Portable and mobile radio frequency communications equipment can affect the UVT Device. The UVT device should not be used in close proximity to or in conjunction with other equipment.

Note:

- Use of devices such as cellular phones, transceivers, or remote control devices that generate radio waves near this product may cause interference with the operation of this product;
- When installing, ensure that this equipment is located as far away from other electrical equipment as possible;
- Use only cables provided or designed by our company;
- Do not attempt to modify this equipment. Modifications to the product may reduce the EMC performance. Modifications include modifications to cables, failure to properly install or position the product, changes to its configuration or components, changes to established operating procedures for the product or accessories, etc.

2. Measures to eliminate EMC problems

If this equipment does cause interference, the user is encouraged to correct the problem by one or more of the following measures:

- Locate other devices as far as possible from this device to reduce electromagnetic interference;
- To reduce electromagnetic interference, you can adjust the position of the equipment and other devices or the angle between them;
- Electromagnetic interference can be reduced by changing the connection point of the UHT device's power cord;
- You can also change the power path to other devices to reduce electromagnetic interference.

Guidance and Manufacturer's Declaration – Electromagnetic Interference

SALUS-FSWT is intended for use in the electromagnetic environment described below. Purchasers or users of SALUS-FSWT must ensure that it is used under such conditions.

Emissions test	Correspondence	Electromagnetic Environment - Guidance
Radio interference GOST CISPR 11-2017	Group 1	SALUS-FSWT uses radio frequency energy only For their internal functions . That's why his radio emission Very low And hardly whether will cause any interference V located near electronic equipment
Radio interference GOST CISPR 11-2017	Class A	SALUS-FSWT suitable For use in everyone institutions except residential , and Maybe be used V household institutions And V those who directly connected To public low voltage networks power supply buildings used For household goals and self-care.
Harmonic radiation GOST IEC 61000-3-2-2021	Class A	
Voltage fluctuations/ flickering radiation GOST IEC 61000-3-3-2015	Compliant	

Guidance and Manufacturer's Declaration - Electromagnetic Emissions

SALUS-FSWT is intended for use in the electromagnetic environment described below.
Purchasers or users of SALUS-FSWT must ensure that it is used under such conditions.

Immunity test IEC 60601	Test Level IEC 60601	Compliance Level	Electromagnetic Environment - Guidance
Electrostatic discharge (ESD) GOST 30804.4.2-2013	± 6 kV contact ± 8 kV atmospheres	± 6 kV contact ± 8 kV atmospheres	Floors should be wood, concrete or ceramic tiles. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transients/pulses GOST IEC 61000-4-4-2016	± 2 kV (For lines power transmission) ± 1 kV (For entrance / exit lines)	± 2 kV (For lines power transmission) ± 1 kV (For entrance / exit lines)	The quality of the power supply should be that of a typical commercial or hospital environment.
Overvoltage GOST IEC 61000-4-5-2017	line ± 1 kV Line to ground ± 2 kV	± 1 kV differential mode ± 2 kV , common mode	The quality of the power supply should be that of a typical commercial or hospital environment.

Voltage dips, short-term interruptions and voltage fluctuations on input power lines. GOST 30804.4.11-2013	About 0.5 cycle <5% UT (>95% reduction in UT) About 5 cycles 40% UT (60% reduction in UT)About 25 cycles 70% UT (30% reduction in UT) About 5 seconds <5% UT (>95% reduction in UT)	About 0.5 cycle <5% UT (>95% reduction in UT) About 5 cycles 40% UT (60% reduction in UT)About 25 cycles 70% UT (30% reduction in UT) About 5 seconds <5% UT (>95% reduction in UT)	Power quality should be that of a typical commercial or hospital environment. If the FSWT user requires continuous operation during power outages, it is recommended that the system be powered by an uninterruptible power supply or battery .
Mains frequency (50/60 Hz) Magnetic field GOST IEC 61000-4-8-2013	30 A/m	30 A/m	Magnetic fields industrial frequencies should be on level characteristic For typical commercial or sick leave environment
NOTE : UT is the power supply voltage			

Guidance and Manufacturer's Declaration - Electromagnetic Immunity

SALUS-FSWT is intended for use in the electromagnetic environment described below. Purchasers or users of SALUS-FSWT must ensure that it is used under such conditions.			
Noise immunity test	Test Level IEC 60601	Compliance Level	Electromagnetic Environment - Guidance
Conductive radio frequency electromagnetic field GOST IEC 61000-4-6-2022	3 V rms 150 kHz – 80 MHz	3 V rms	Portable or mobile communications equipment should not be used on any cable section closer than this separation distance calculated

Resistance to RF fields GOST IEC 61000-4-3- 2016	3 V/m 80 MHz – 2.5 GHz	3 V/m	from the equation applicable to the transmitter frequency. 80 MHz – 800 MHz $d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz – 2.5 GHz Where P is the maximum rated output power of the transmitter in watts (W) as specified by the transmitter manufacturer, and d is the recommended distance in meters (m). The field strength of a fixed RF transmitter, as determined from an electromagnetic field survey, must be less than the compliance level for each frequency range. B Interference may occur in the vicinity of equipment marked with the following symbol: 
NOTE 1: At 80 MHz and 800 MHz applies more high frequency range .			
NOTE 2: These recommendations can apply Not in everyone situations On spreading electromagnetic radiation influences absorption And reflection from structures , objects And people			
^a Field strengths from fixed transmitters such as radio (cellular/cordless) telephones and base stations for land mobile radios, amateur radio, AM and FM radio and television broadcasts cannot be theoretically predicted with accuracy. To assess the electromagnetic environment created by a stationary radio transmitter, an electromagnetic field survey should be considered. If the measured field strength at the location where SALUS-FSWT is used exceeds applicable RF compliance level, monitoring is required to ensure that the SALUS-FSWT is operating properly. If abnormalities are observed, additional measures such as orientation or adjustment of the SALUS-FSWT may be required. b. The field strength must be less than [V1] V/m in the frequency range from 150 kHz to 80 MHz.			

Recommended distance between portable and mobile radio communication devices and SALUS-FSWT

SALUS-FSWT is intended for use in an electromagnetic environment in which radiated radio frequency interference is controlled. The purchaser or user of the SALUS-FSWT can help prevent electromagnetic interference by maintaining a minimum separation distance between SALUS-FSWT portable and mobile RF communications equipment (transmitters) as recommended below, consistent with the maximum communications output power.			
Rated maximum output power of the transmitter, W	Separation distance depending on transmitter frequency		
	150 kHz - 80 MHz	80 MHz - 800 MHz	800MHz - 2.5GHz
	$d = \left[\frac{3.5}{V_1} \right] \sqrt{P}$ Where $V_1 = 3$	$d = \left[\frac{3.5}{E_1} \right] \sqrt{P}$ Where $E_1 = 3$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ Where $E_1 = 3$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

The recommended distance d (m) for a transmitter rated for a maximum output power not specified above can be determined using the formula applicable to the frequency of the transmitter, where P is the power in watts (W) is the maximum rated output power of the transmitter

NOTE 1: Apply higher frequency range 80 MHz and 800 MHz.

NOTE 2: These guidelines may not apply in all situations. The propagation of electromagnetic radiation is affected by absorption and reflection from structures, objects and people.

The manufacturer assumes no responsibility for interference caused by the use of a non-recommended connecting cable or by unauthorized modification or alteration of the device. Unauthorized modifications or alterations to the UVT equipment may void the user's authority to operate the UVT equipment.

List of regulatory documentation used by the manufacturer to confirm compliance:

EN 1041 , Information provided by the manufacturer of medical devices
ISO 13485 , Medical devices – Quality management systems – Requirements for regulatory purposes
ISO 14971 , Medical devices – Application of risk management to medical devices
ISO 15223-1 , Medical devices - Symbols to be used on medical device labels, markings and information to be provided - Part 1: General requirements
IEC 60601-1 , Medical electrical equipment – Part 1: General requirements for basic safety and essential performance
IEC 60601-1-2 , Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance - Related standards: Electromagnetic compatibility: Requirements and tests
IEC 60601-1-6 , Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Associated standard: Usability
IEC 60601-2-10 , Medical electrical equipment – Part 2-10: Particular requirements for the safety of nerve and muscle stimulators
IEC 62304 , Medical device software – Software life cycle processes
IEC 62366 , Medical devices – Application of serviceability
MEDDEV 2.7/1 Revision 4 , CLINICAL EVALUATION: GUIDANCE FOR MANUFACTURERS AND NOTIFIED BODIES IN ACCORDANCE WITH DIRECTIVES 93/42/ EEC and 90/385/ EEC
ISO 10993-1 , Biological evaluation of medical devices – Part 1: Evaluation and testing as part of the risk management process
ISO 10993-5 , Biological evaluation of medical devices – Part 5: Cytotoxicity tests by in vitro
ISO 10993-10 , Biological evaluation of medical devices – Part 10: Tests for skin irritation and sensitization
EN 980 - Symbols used on the labeling of medical devices

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