POLYMAILLE®

- FR** Prothèses vasculaires en polyester tricoté
- EN** Knitted-polyester vascular prostheses
- DE** Gefäßprothesen aus gestricktem Polyester
- ES** Prótesis vasculares de poliéster tejido
- IT** Protesi vascolari in poliestere intrecciato
- PT** Próteses vasculares em poliéster tricotado
- NL** Vaatprothesen uit gebreide polyester
- PL** Protezy naczyńiowe z poliestru trykotowego
- DA** Vaskulære proteser af strikket polyester
- SV** Av polyester stickade kärlproteser
- NO** Vaskulær protese i strikket polyester
- FI** Polyesteritrikoiset verisuoniproteesit
- EL** Αγγειακά μοσχεύματα από πλεκτό πολυεστέρα
- TR** Örülü polyester vasküler protezler
- RO** Proteză vasculară din poliester tricotat
- HU** Kötött poliészter érprotézisek
- RU** Протезы сосудов из полиэфирного трикотажа
- CS** Cévní protézy z pleteného polyesteru
- AR** الوصلات الوعائية من البوليستر المحبوك
- BG** Съдови протези от тъкан полиестер

PEROUSE

STERILE R FR : stérilisé par irradiation EN : sterilized by irradiation DE : sterilisiert durch Bestrahlung ES : esterilizado por irradiación IT : sterilizzato con irradiazione PT : esterilise pela irradiação NL : gesteriliseerd per uitstraling PL : sterylizacji metodą napromieniowania DA : steriliseret ved bestråling SV : steriliserats med strålning NO : sterilisert ved bestråling FI : teriloitu säteilyty EL : αποστειρωμένος από την ακτινοβολία TR : ışınlama ile steril RU : стерилизир при ирадиере HU : sterilizáláni besugárzás RO : prosterilizovano obлъчением CS : sterilizovaný zářením AR : تعقيمها من تشعيع BG : стерилизирани чрез облъчване	 FR : Lire la notice EN : See instruction for use DE : Anleitung lesen ES : Consulte las instrucciones de uso IT : Leggere attentamente le istruzioni per l'uso PT : Consulte as instruções de utilização NL : Lees de handleiding PL : Zobacz instrukcję obsługi DA : Se brugsanvisning SV : Se bruksanvisningen NO : Les bruksanvisningen FI : Lue käyttöohje EL : Βλ. οδηγίες χρήσης TR : Kullanma kılavuzuna bakınız RU : См. инструкцию по использованию UK : Дивіться інструкцію з використання CS : Viz návod k použití RO : Consultați instrucțiunile de utilizare HU : Lásd a használati útmutatót AR : الإطلاع على نشرة التعليمات BG : Вижте указанията
 FR : Usage unique, EN : Single use, DE : Zum einmaligen Gebrauch ES : De un solo uso IT : Monouso PT : Utilização única NL : Eenmalig gebruik PL : Do jednorazowego użytku DA : Engangsbrug SV : Engångsbruk NO : Engangsbruk FI : Kertäkäyttö EL : Μίας χρήσης TR : Tek kullanımlık RU : Одноразового применения UK : Одноразового використання CS : Na jedno použití RO : De unică folosință HU : Egyszeri használatra AR : تستخدم مرة واحدة BG : За еднократна употреба	 FR : Conserver à l'abri de la lumière du soleil EN : Keep away from sunlight DE : Unter Lichtabschluss aufbewahren ES : Conservar en un lugar protegido de la luz solar IT : Conservare al riparo dalla luce del sole PT : Manter abrigado da luz solar NL : Beschermen tegen het zonlicht PL : Przechowywać z daleka od słońca DA : Beskyttes mod sollys SV : Skyddas från solljus NO : Lagres beskyttet fra sollys FI : Säilyttävä suojassa auringonvalolta EL : Να φυλάσσεται μακριά από το ηλιακό φως TR : Güneş ışığından uzak tutunuz RU : Хранить в месте, защищенном от солнечных лучей UK : Зберігати в місці, яке захищено від сонячних променів CS : Zabránit přístupu slunečního světla RO : A se păstra la adăpost de lumina soarelui HU : Napfénytől távol tartandó AR : يحفظ بعيداً عن أشعة الشمس BG : Да се пазят от пряка слънчева светлина
 FR : Ne pas utiliser si l'emballage est endommagé EN : Do not use if the packaging is damaged DE : Nicht verwenden, wenn die Verpackung beschädigt ist ES : No utilizar si el embalaje está deteriorado IT : Non utilizzare se la confezione è danneggiata PT : Não utilizar se a embalagem se encontrar danificada NL : Niet gebruiken als de verpakking beschadigd is PL : Nie używać jeśli opakowanie jest uszkodzone DA : Må ikke anvendes, hvis emballagen er beskadiget SV : Får inte användas om förpackningen är skadad NO : Skal ikke brukes hvis emballasjen er skadet FI : Ei saa käyttää, jos pakkaus on vioittunut EL : Να μην χρησιμοποιείται αν η συσκευασία έχει υποστεί ζημιά TR : Ambalaj hasar görmüşse kullanılmayınız RU : Не использовать в случае повреждения упаковки UK : Не використовувати у разі пошкодження упаковки CS : Nepoužívejte, pokud je porušen obal RO : A nu se utiliza dacă ambalajul este deteriorat HU : Ne használja fel, ha a csomagolás sérült AR : لا تستخدم في حالة وجود أي تلف بالعبوة BG : Да не се употребява, ако е нарушена целостта на	 FR : Craint humidité EN : Protect from humidity DE : Vor Feuchtigkeit schützen ES : Proteger de la humedad IT : Teme l'umidità PT : Manter a humidade NL : Kan niet tegen vocht PL : Trzymać z daleka od wilgoci DA : Beskyttes mod fugt SV : Skyddas från fukt NO : Beskyttes mot fuktighet FI : Suojattava kosteudelta EL : Να προστατεύεται από την υγρασία TR : Nemi konulunuz RU : Поддерживайте влажность UK : Підтримувати вологість CS : Udržujte vlhkost RO : A se proteja de umezeală HU : Száras helyen tartandó AR : يحفظ بعيداً عن الرطوبة BG : Да се пазят от влага

LOT FR : Numéro de lot EN : Batch code DE : Chargennummer ES : Número de lote IT : Numero di lotto PT : Número de lote NL : Lotnummer PL : Numer partii produktu DA : Batch-kode SV : Lotnummer NO : Lotnummer FI : Eränumero EL : Κωδικός παρτίδας TR : Parti kodu RU : Номер партии UK : Номер партії CS : Číslo šarže RO : Număr lot HU : Adagszám AR : رقم الشحنة BG : Партида номер	SN FR : Numéro de série EN : Serial number DE : Seriennummer ES : Número de serie IT : Numero di serie PT : Número de série NL : Serienummer PL : Numer serii DA : Serienummer SV : Serienummer NO : Serienummer FI : Sarjanumero EL : Αριθμός σειράς TR : Seri numarası RU : Номер серии UK : Номер серії CS : Sériové číslo RO : Număr serie HU : Sorozatszám AR : رقم التسلسل BG : Серийен номер
 FR : À utiliser avant la date EN : Use by DE : Verfallsdatum ES : Utilizar antes de IT : Utilizzare entro PT : Utilizar antes de NL : Te gebruiken vóór PL : Do użycia przed DA : Anvendes inden SV : Används före NO : Brukes før FI : Käytettävä ennen EL : Χρηστέος έως TR : Son kullanma tarihi RU : Использовать до UK : Використовувати до CS : Použítelné do RO : A se utiliza înainte de HU : Minőségét megőrzi AR : تاريخ انتهاء الصلاحية BG : Да се използва преди	REF FR : Numéro de référence EN : Reference number DE : Bestellnummer ES : Número de referencia IT : Numero di catalogo PT : Referência de catálogo NL : Catalogusnummer PL : Numer katalogowy DA : Referencenummer SV : Referensnummer NO : Henvisningsnummer FI : Viitenumero EL : Αριθμός αναφοράς TR : Referans numarası RU : Номер по каталогу UK : Номер за каталогом CS : Referenční číslo RO : Număr de referință HU : Hivatkozási szám AR : كود المنتج BG : Каталоген номер
 FR : Emballage recyclable EN : Recyclable packaging DE : Verpackung wiederverwertbar ES : Embalaje reciclable IT : Confezione riciclabile PT : Embalagem reciclável NL : Recycleerbare verpakking PL : Opakowanie do użytku wtórego DA : Emballagen kan genanvendes SV : Förpackningen kan återvinas NO : Resirkulerbar emballasje FI : Kierrätettävä pakkaus EL : Ανακυκλώσιμη συσκευασία TR : Geri dönüştürülebilir ambalaj RU : Упаковка подлежит вторичной переработке UK : Упаковка підлягає вторинному переробленню CS : Recyklovatelný obal RO : Ambalaj reciclabil HU : Újrahasznosítható csomagolás AR : عبوة قابلة لإعادة التدوير BG : Опаковка за рециклиране	

PEROUSE MEDICAL

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INSTRUCTIONS AND INDICATIONS

I/ DESCRIPTION

POLYMAILLE® vascular prostheses are constructed from knitted polyester. The following models are currently available:

- POLYMAILLE® C and POLYMAILLE® Extra Thin, with or without SHR® radial reinforcement, impregnated with collagen of bovine origin;
- POLYMAILLE® N, non-impregnated.

The POLYMAILLE® N vascular prosthesis requires pre-clotting prior to insertion as it has a water permeability profile of 1000 ml/cm²/min at the pressure of 120 mmHg.

Collagen impregnation of the POLYMAILLE® C and POLYMAILLE® Extra Thin vascular prostheses ensures water permeability profile of less than 5 ml/cm²/min at 120 mmHg and allows direct implantation without the need for pre-clotting, thereby reducing operating time and blood loss.

The POLYMAILLE® Extra Thin prosthesis has a wall thickness of about 0.32 mm. In addition, its compliance approaches that of a native artery, which makes it pulsatile.

The non-divisible radial reinforcement of POLYMAILLE® prosthesis with SHR® (Stretch Helicoidal Reinforcement) technology ensures improved protection from external compressive stresses. This helical reinforcement matches the corrugated profile the prosthesis, providing a "stretch" effect making the reinforced portion "stretchable", and reducing the anastomotic stresses, while ensuring improved kink resistance.

II/ INDICATIONS

POLYMAILLE® vascular prostheses are indicated for replacement or bypass of arteries presenting aneurysm or obliterative arterial disease.

Their indication is restricted to abdominal and peripheral surgery not involving the crossing of the knee joint.

III/ CONTRAINDICATIONS

POLYMAILLE® vascular prostheses are not indicated for surgery of the thoracic aorta, coronary bypass or arteriovenous access.

The POLYMAILLE® C and POLYMAILLE® Extra Thin prostheses are contraindicated in patients with known hypersensitivity to collagen. The POLYMAILLE® N prosthesis is contraindicated in patients on high-dose heparin and patients with clotting disorders.

IV/ POSSIBLE COMPLICATIONS

The most commonly encountered complications following implantation of polyester vascular prostheses include, but are not restricted to: infection, periprosthetic seroma, thrombosis, pseudoaneurysms, periprosthetic haematoma.

V/ WARNINGS AND PRECAUTIONS

POLYMAILLE® vascular prostheses are for single use only. DO NOT REUSE.

POLYMAILLE® vascular prostheses are sterilised by irradiation. DO NOT RESTERILIZE.

DO NOT USE IF THE PACKAGE IS DAMAGED

INSTRUCTIONS AND INDICATIONS

DO NOT USE POLYMAILLE® VASCULAR PROSTHESES AFTER THEIR EXPIRY DATE.

Store protected from light, cold, heat (10°C < T < 35°C) and humidity (RH < 80%).

DO NOT PRE-CLOT THE POLYMAILLE® C and POLYMAILLE® Extra Thin VASCULAR PROSTHESES: they are impregnated with collagen.

PRE-CLOTting IS ESSENTIAL BEFORE INSERTION OF THE POLYMAILLE® N PROSTHESES.

The SHR® RADIAL REINFORCEMENT IS NON-DIVISIBLE, DO NOT REMOVE IT.

For prostheses with SHR® reinforcement, removing a portion of the non-divisible reinforcement poses significant risks of leakage, weakening and failure of the prosthetic wall.

For prostheses with SHR® reinforcement, take into account the limited elongation of the reinforced portion, which is expected to occur in the medium term because of the flattening of the corrugated peaks of the initial corrugation.

Prostheses contaminated with blood and/or other body fluids prior to surgery must not be used.

Atraumatic clamps should be used or the jaws of the clamps should be enclosed in soft protective covers. Tighten as little as possible.

Avoid excessive tension on the prosthesis during insertion, including for prostheses with SHR® reinforcement.

Avoid excessive length and axial torsions (twisting), which may result in kinking.

Suturing should be performed using only tapered needles; sharp-edged needles may damage the prosthesis.

For prostheses with SHR® reinforcement, perform anastomoses only on non-reinforced portions.

VI/ INSTRUCTIONS FOR USE

Vascular prostheses should only be inserted by surgeons trained in vascular surgery techniques.

On removal from their protective packaging, the prostheses are as follows:

- POLYMAILLE® N comes packed in a TYVEK® peel-open sachet containing a sterile blister pack holding the prosthesis.

Carefully peel open the TYVEK® sachet and hand the contents only to theatre staff wearing sterile clothing. The second pack (blister) is sterile and should be opened by the surgical team. New sterile gloves should be used to handle the prosthesis.

- POLYMAILLE® C and POLYMAILLE® Extra Thin come in a peel-open aluminium sachet containing two packs (a peel-open Tyvek® sachet containing a sterile blister pack holding the prosthesis).

ATTENTION: DO NOT OPEN THE ALUMINIUM SACHET ONTO THE STERILE DRAPE.

This sachet is simply a protective package and its contents are not sterile.

Peel open the aluminium sachet and proceed as with POLYMAILLE® N.

The diameter of the prosthesis must correspond to the diameter of the patient's vascular lumen. Allowance must be made for the patient's weight, the operating position and the degree of movement likely to occur at the implantation site when determining the correct

INSTRUCTIONS AND INDICATIONS

length of the prosthesis to be inserted.

VII/ TRACEABILITY

Self-adhesive labels and a patient card are provided in each packaging to allow the required follow-up of implanted prostheses. One label must be pasted onto the patient's medical record and another onto the patient card.

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INSTRUCTIONS AND INDICATIONS

I/ DESCRIPTION

POLYTHESE® vascular prostheses are woven polyester vascular grafts. Three models are available:

- POLYTHESE® IC and ICT, impregnated with collagen of bovine origin;
- POLYTHESE® IC-3GL, impregnated with collagen of bovine origin;

Collagen impregnation of POLYTHESE® IC, IC-3GL and ICT vascular prostheses ensures a water permeability of less than 10 ml/cm²/min at 120 mmHg and allows direct implantation of the prosthesis without the need for pre-clotting, thereby reducing operating time and blood loss.

POLYTHESE® IC-3GL prostheses have 3 equidistant longitudinal guide lines.
FORMALDEHYDE/GLUTARALDEHYDE-FREE

II/ INDICATIONS

POLYTHESE® vascular prostheses are indicated for replacement or bypass of arteries presenting aneurysm or arterial disease.

Their indication is restricted to thoracic (POLYTHESE® ICT), abdominal and peripheral surgery not involving the crossing of the knee joint.

The POLYTHESE® IC-3GL prostheses are especially indicated in aortic root and ascending thoracic aorta surgery.

III/ CONTRAINDICATIONS

POLYTHESE® vascular prostheses are not indicated for coronary bypass or arteriovenous access.

POLYTHESE® IC, IC-3GL and ICT prostheses are contraindicated for patients with known hypersensitivity to collagen.

IV/ POSSIBLE COMPLICATIONS

The most commonly encountered complications following implantation of polyester vascular prostheses include, but are not restricted to: infection, periprosthetic fluid collection, thrombosis, pseudoaneurysms, periprosthetic haematoma, leaking, rupture, lymphatic complications, venous thrombosis, pulmonary embolism.

V/ WARNINGS AND PRECAUTIONS

The afore-mentioned warning messages are to be strictly respected; for the patient, their non-compliance can result in significant blood loss, disability of a limb or even death.

- **POLYTHESE® vascular prostheses are for single use only. DO NOT REUSE.** POLYTHESE® prostheses must never be reimplanted. As they are contaminated with blood, a used vascular prosthesis constitutes a biological risk (infectious risk). In addition, PEROUSE MEDICAL cannot guarantee proper functioning of the device after explantation (water-tightness of the prosthesis; integrity of the collagen coating; weakening of the textile structure, particularly at the areas of anastomosis; level of burst strength; level of suture retention strength, etc.)
- **POLYTHESE® vascular prostheses are sterilised by irradiation. DO NOT RESTERILIZE.**
- **DO NOT USE THE IMPLANT IF THE PACKAGE IS DAMAGED.**
- **DO NOT USE POLYTHESE® VASCULAR PROSTHESES AFTER THEIR EXPIRY DATE.**

INSTRUCTIONS AND INDICATIONS

- Store away from light, cold, heat (10°C < T < 35°C) and humidity (RH < 80%).
- **DO NOT PRE-CLOT POLYTHESE® IC, IC-3GL AND ICT PROSTHESES: they are impregnated with collagen.**
- Prostheses contaminated with blood and/or other body fluids prior to surgery must not be used.
- Do not use scissors to cut the prosthesis to the desired length.
- The use of a thermocutter can prevent fraying of the extremities of the prosthesis.
- Atraumatic clamps should be used or the jaws of the clamps should be enclosed in soft protective covers. Tighten as little as possible.
- Avoid excessive tension on the prosthesis during implantation.
- Avoid excessive length and axial torsions (twisting), which may result in kinking.
- Suturing should be performed using only tapered needles; sharp-edged needles may damage the prosthesis.

VI/ INSTRUCTIONS FOR USE

- Vascular prostheses should only be implanted by surgeons trained in vascular surgery techniques.
- On removal from their outer packaging, the prostheses are as follows: peel-open aluminium pouch containing the double packaging (external blister pack containing a second sterile blister pack holding the prosthesis.).

WARNING: DO NOT OPEN THE ALUMINIUM POUCH ON THE STERILE DRAPE.

It is simply a protective package and its contents are not sterile.

Peel open the aluminium pouch then carefully open the first blister back and hand the contents only to people wearing sterile clothing. The second blister pack is sterile and should be opened by the surgical team. New sterile gloves should be used to handle the prosthesis.

- The diameter of the prosthesis must correspond to the diameter of the patient's vascular lumen. The patient's weight must be considered the operating position and the degree of movement likely to occur at the implantation site when determining the correct length of the prosthesis to be implanted.

VII/ TRACEABILITY

Self-adhesive labels and a patient card are provided in each packaging to allow the required follow-up of implanted prostheses. One label must be stuck onto the patient's medical record and another onto the patient card.

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[Signature]

POLYMAILLE® FLOW PLUS HEPARIN

Heparin-coated ePTFE vascular prosthesis

Product description

POLYMAILLE® FLOW PLUS HEPARIN is a vascular prosthesis made of expanded Polytetrafluorethylene (ePTFE) and featuring a heparin-coated surface. The prosthesis comprises a single-layer prosthesis corpus with a thin ePTFE mantle wound around it on the outside. The lumen surface of the POLYMAILLE® FLOW PLUS HEPARIN prosthesis is coated with highmolecular heparin (pig mucosa).

This vascular prosthesis is available in the two wall thicknesses: „thin wall“ and „standard wall.“ In areas of potentially high levels of compression and bending pressure (interarticular and extraanatomical reconstruction), it is also available fitted with a removable external spiral reinforcement made of Polytetrafluorethylene (PTFE).

Four patient stickers are enclosed in the package for purposes of patient documentation.

Indications

The POLYMAILLE® FLOW PLUS HEPARIN is indicated in arterial reconstruction, primarily in peripheral blood vessels. It can also be used extraanatomically in femoro-femoral and axillofemoral vascular reconstruction.

ООО ЦИРМИ  INFO@CIRMI.RU  WWW.GOSZDRAVNADZOR.RU

The standard wall POLYMAILLE® FLOW PLUS HEPARIN prosthesis is also indicated as an arterio-venous shunt prosthesis in haemodialysis.

Contraindications

The clinical data on record to date do not provide a sufficient basis for conclusions concerning the use of ePTFE vascular prostheses in the following indications:

Reconstruction

- of the pulmonary arteries
- of the aorta ascendens
- of the coronary arteries
- of the carotid arteries (a. carotis communis, externa and interna)
- of the cerebral arteries
- of the arteria brachialis and cephalica
- of the cardiac and pulmonary veins as well as
- of the vena cava inferior and/or superior

The POLYMAILLE® FLOW PLUS HEPARIN should also not be implanted in patients

- who show an immunological reaction to ePTFE or
- who show a hypersensitive reaction to heparin (heparin allergy), including patients with HIT type II in their medical histories.

Interactions with other medical drugs

Administration of protamine influences the biocompatibility properties of heparinized surfaces. Administration of protamine should therefore not be carried out in the immediate vicinity of the coated prosthesis. Minimization of protamine treatments and staged dosing of protamine is recommended. Interactions with other medical drugs are unknown to date.

Information and precautions

Please bear in mind that serious complications may arise to the patient if following notes are not observed.

1. The POLYMAILLE® FLOW PLUS HEPARIN ePTFE vascular prosthesis should only be implanted by persons trained in the surgical implantation technique required for the material ePTFE.
2. Anticoagulation measures should not be reduced more than usual.
3. Preoperative diagnostics are recommended to detect any increase in risk of haemorrhaging in patients receiving oral anticoagulant therapy (e.g. acetylsalicylic acid, warfarin).
4. Do not use this product after the expiry date. Please consult your supplier in such cases.
5. The POLYMAILLE® FLOW PLUS HEPARIN ePTFE prosthesis must not be resterilized.
6. The prosthesis should not be exposed to high levels radiation, due to which the ePTFE apart.
7. Do not allow the prosthesis to come into contact with reagents containing oil, alcohol or aqueous solutions. Also, do not attempt to preclot the prosthesis.

8. Avoid excessive processing or manipulating of the prosthesis, since otherwise blood, serum or plasma could be pressed through the prosthesis wall.
9. Use atraumatic clamps with soft, cushioned ends.
10. Always use non-absorbable suture material. The suture seam on vessels should not measure less than 2 mm.
11. Thin wall ePTFE vascular prostheses cannot be recommended for extraanatomical, axillo-femoral vascular prosthetics, since there is a risk of anastomosis rupture under heavier stress loads.
12. POLYMAILLE® FLOW PLUS HEPARIN prostheses with a spiral reinforcement over their entire length are not indicated for use as arterio-venous shunt prostheses. Do not puncture in the area of reinforcement. This would involve an increased risk of embolism and/or formation of pseudoaneurysms.
13. Do not use the prosthesis or parts thereof as a patch. The product is not designed so as to be suitable for indications other than those listed or to achieve good results in such indications. On the contrary, such use may lead to adverse events.
14. POLYMAILLE® FLOW PLUS HEPARIN is a single-use product that has been sterilized with EtO. It may not be resterilized or reused. Risks arising from resterilization and reuse include infection, loss of function, hemolysis, allergic reaction in the patient. The product must not be used, if the package is damaged or if the package has been opened outside of the sterile area.

- It is recommended that the same width of anastomosis prosthesis / tissue suture seam be maintained on both sides. This seam margin should not be less than 2 mm in reconstruction of vessels.
- Avoid excessive tension on the anastomosis, but also do not leave excessive suture gaps.
- With the axillo-femoral bypass, the recommended position for the proximal anastomosis is as far to medial as possible within the first third of the a. axillaris. The anastomosis should be performed with the patient's arm extended.
- Minimize the risk of puncture channel bleeding by following the curve of the needle, holding the suture thread at a 90° angle as far as possible and not pulling too hard on the suture.

Thrombectomy:

- POLYMAILLE® FLOW PLUS HEPARIN ePTFE vascular prostheses are not radially elastic. When balloon catheters and/or angioplasty catheters are to be used postoperatively in the prosthesis lumen, make absolutely sure the outside diameter of the inflated balloon corresponds to the inside diameter of the prosthesis. Overinflation and/or use of the wrong size may result in dilatation of and/or damage to the prosthesis.

POLYMAILLE® FLOW PLUS HEPARIN ePTFE as an a.v. shunt prosthesis in haemodialysis:

- To reduce haematoma formation, we recommend not puncturing the shunt until about 2 weeks after the implantation.

- Make sure the dialysis cannula is first inserted into the prosthesis at a 45° angle, then guided parallel to the prosthesis.
- Do not puncture the same locus or area repeatedly, since this may lead to formation of haematomas and pseudo neurysms.

Storage and packaging

The POLYMAILLE® FLOW PLUS HEPARIN ePTFE vascular prosthesis is an EtO-sterilized, disposable (single-use) product that is supplied in a sterile package and will remain sterile as long as the package has not been opened or damaged.

Store in a cool, clean, dry environment at temperatures between 10°C and 40°C.

Exchange guarantee

If the package has been accidentally opened or damaged, or if damage to the package or product has occurred during transport, send the complete product in all of the packaging material back to your supplier without delay. The supplier will quickly exchange the product accordingly.

GUARANTEE AND WARRANTY

PEROUSE MEDICAL GmbH guarantees that this product complies with the applicable legal requirements, in particular as stipulated by the EC Directive for Medical Products, 93/42/EEC. PEROUSE MEDICAL GmbH will therefore, should the product be shown to exhibit material or processing defects before the expiry date, which defects resulted during production or packaging,

- replace the product with one of identical construction or
- with a product equivalent to the original as to function.

All of the following conditions must be met prior to avilment of the guarantee:

- The product must not have been used after its expiry date.

- b) First manifestation of the defect must not have been after the product expiry date.
- c) Use of the product revealing the defect must have been the first use of a product licensed for single use only or intended for single use only by PEROUSE MEDICAL GmbH.
- d) The defect must have been reported to PEROUSE MEDICAL GmbH or its representative without delay.
- e) The defective product must have been sent back to PEROUSE MEDICAL GmbH or its representative immediately after the defect was reported.

The above guarantee shall not affect or change any of the required stipulations deriving from the seller's warranty and product liability legislation and shall apply in addition to them.

Otherwise, the guarantee applies only within the scope defined above. Further claims may not be made based on the guarantee. This guarantee, as well as the above-mentioned legal warranty and liability claims, shall apply in place of and to the exclusion of any other commitments or assurances made by PEROUSE MEDICAL GmbH, whether oral or written, express or tacit in nature. None of our representatives or employees has the right to make any commitments or assurances to the contrary or to enter into agreements to the contrary.

In view of the many factors beyond the control of PEROUSE MEDICAL GmbH such as shipping, storage and handling by the user, improper use of or non-compliance with the instructions for use or warnings, PEROUSE MEDICAL GmbH is in particular not liable for damages of any kind, resulting either directly or indirectly from the product or its use, unless such damages are proved to have been caused by a defect or defects in the product.

In particular, PEROUSE MEDICAL GmbH assumes no liability for damages resulting from resterilization or reuse of products approved for single use only, or which PEROUSE MEDICAL GmbH designates as being for single use.

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ООО ЦРПМ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

[Handwritten signature]

Quality and Regulatory Affairs Manager

English	P. 4
Français	P. 12
Deutsch.....	P. 21
Español	P. 30
Italiano.....	P. 38
Português	P. 47
Polski	P. 56
Türkçe	P. 65

POLYMAILLE® FLOW

ePTFE Vascular Prosthesis

Product description

POLYMAILLE® FLOW is a vascular prosthetic graft made of expanded polytetrafluoroethylene (ePTFE). A thin ePTFE wrap is wound around the single-layer body of the graft.

The graft is available in two wall thicknesses, „thin wall“ and „standard wall“. It is also available with a removable external helical Polytetrafluoroethylene (PTFE) reinforcement for sections in which the risk of compression and kinking is high (reconstructions that traverse a joint or extraanatomical reconstructions).

Four stickers for patient documentation are enclosed in the packaging.

Indications

The POLYMAILLE® FLOW is indicated in cases of arterial vasoreconstruction, primarily in the peripheral vascular region. It can also be used in extraanatomical reconstructions – femorofemoral and axillofemoral.

Standardwall POLYMAILLE® FLOW grafts are also indicated for use as arteriovenous shunt prostheses in haemodialysis.

Contraindications

The available clinical data are insufficient to support definitive conclusions concerning the use of ePTFE vascular grafts in the following areas of application:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Reconstruction

- of the pulmonary arteries
- of the aorta ascendens
- of the coronary arteries
- of the carotid arteries (A. carotis communis, externa and interna)
- of the cerebral arteries
- of the arteria brachialis and cephalica
- of the cardiac and pulmonary veins and
- of the vena cava inferior and/or superior

The POLYMAILLE® FLOW should also not be implanted in patients who show an immunological hyperreaction to ePTFE.

Precautions and warnings

Failure to comply with these precautions and warnings may result in serious medical complications.

- The POLYMAILLE® FLOW ePTFE vascular prosthetic graft should only be implanted by persons with special training in the surgical implantation of the material ePTFE.
- Do not use this product after the expiry date. Please ask your supplier in such cases.
- The graft must not be exposed to strong radiation. ePTFE decomposes when exposed to high doses of radiation.
- Do not allow the graft to come into contact with reagents containing oil, alcohol or aqueous solutions. Do not attempt to preclot the graft.

- Keep processing and manipulation of the graft to a minimum, since otherwise blood, serum or plasma could be pressed through its wall.
- Use atraumatic clamps with soft, cushioned ends.
- Always use non-absorbable suture material. The suture seam should not measure less than 2 mm on vessels.
- Thin-walled ePTFE vascular prostheses are generally not recommended for extraanatomical, axillofemoral vessel prosthesis, since the anastomosis may rupture under heavy stress loads.
- POLYMAILLE® FLOW grafts with helical reinforcement over their entire length are not indicated for use as arteriovenous shunt prostheses. In the graft models with central helical reinforcement, blood access is not permissible in the area of the reinforcement. There is an increased risk of embolism and/or pseudoaneurysms in both cases.
- Do not use this prosthetic graft or parts of it as a patch. The product design is not suitable or beneficial for applications other than those indicated and undesirable events may result from such misuse.
- POLYMAILLE® FLOW is a single-use product that has been sterilized with EtO. It may not be resterilized or reused. Risks arising from resterilization and reuse include infection, loss of function, hemolysis, allergic reaction in the patient. The product must not be used, if the package is damaged or if the package has been opened outside of the sterile area.

Adverse reactions

As with any surgical implant, prosthetic infection necessitating explantation is a possibility. Maintenance of stringent aseptic surgical conditions is therefore of the utmost importance.

Other potential complications include thrombosis, suture haemorrhaging along the anastomosis, suture rupture, elongation, seroma formation, development of pseudoaneurysms and haematomas.

Risks involved in use as an arteriovenous shunt graft: steal syndrome, skin erosion.

Implantation procedure

Artificial vascular prostheses must always be implanted by surgeons trained in and familiar with the specific techniques of such operations. Please also observe the following points when implanting the POLYMAILLE® FLOW ePTFE vascular prosthetic graft:

Preparation:

- When opening the packaging, please be aware that it is a double-blister. The inside blister only may be brought into the sterile field.
- Touch the graft with sterile, unpowdered gloves only. Use no sharp or pointed instruments and atraumatic vascular clamps only.
- When selecting the diameter, make sure the dimension matches the lumen of the original physiological blood vessel as closely as possible.

- Before cutting the prosthesis to size, first remove a suitable length of the spiral reinforcement along the outside from each end as follows: Take hold of the end of the spiral with the fingers or an atraumatic instrument, then slowly and carefully unwind the spiral reinforcement, holding it at a 90° angle to the prosthesis.
- Cut the graft to size use a sharp scissor or scalpel. Since ePTFE is not elastic in the longitudinal direction, make sure the cut length is correct, taking into account the patient's position, body weight and ensuring complete freedom of movement of the extremity or extremities. If the graft is cut off too short, the suture holes may be stretched, leading to suture haemorrhaging. Cutting it too long may result in kinks.

Placement using tunneler:

- Helix-reinforced grafts require a tunneler with a bullet tip 1 - 2 mm larger than the inside diameter of the graft to facilitate unobstructed passage and placement. Before pulling the graft through the tunnel, replace the bullet tip with one that corresponds exactly to the inside diameter of the graft.
- When using grafts without external reinforcement, the diameter of the tunneler should be the same as the inside diameter of the graft to be implanted.

Anastomosis:

- Use only atraumatic clamps with soft, cushioned ends.
- For suturing use a non-cutting round needle of suitable size and a corresponding non-absorbable monofilament

- The suture seam along the anastomotic graft tissue should have the same width on both sides. This seam margin should not be less than 2 mm in reconstructions of vessels.
- Avoid excessive tension on the anastomosis as well as inappropriate suture gaps.
- For axillofemoral bypasses, the recommended location for the proximal anastomosis is as far medial as possible in the first third of the a. axillaris. Perform the anastomosis procedure with the patient's arm outstretched.
- Minimize the risk of suture haemorrhaging by following the curve of the needle, holding the suture thread at a 90° angle as far as possible and never pulling hard on the suture seam.

Thrombectomy:

- POLYMAILLE® FLOW ePTFE vascular prosthetic grafts are not radially elastic. For postoperative placement of balloon catheters and/or angioplasty catheters in the graft lumen you must therefore make sure the outside diameter of the inflated balloon is the same as the inside diameter of the graft. Excessive inflation and/or use of the wrong size may result in dilatation of and/or damage to the prosthesis.

POLYMAILLE® FLOW ePTFE vascular prosthetic graft used as an arteriovenous shunt graft in haemodialysis:

- To reduce haematoma formation, blood access (catheterization) through the shunt should not be carried out until about 2 weeks after the implantation procedure.

- Make sure the dialysis cannula is inserted at a 45° angle to the graft, then position it parallel to the prosthesis.
- Do not cannulate the same insertion point twice, since repeated puncturing in the same area may result in formation of haematomas and pseudoaneurysms.

Storage and packaging

The POLYMAILLE® FLOW ePTFE vascular prosthetic graft is an EO-sterilized, single-use product that is supplied in a sterile package. It will remain sterile as long as the package is neither opened nor damaged.

Store in a cool, clean, dry place.

Replacement warranty

Please ship the product with all packaging and accessories back to your supplier without delay in case of accidental opening of or damage to the packaging or if the packaging, or indeed the product, has been damaged in transit. The supplier will see to it that the items are replaced immediately.

WARRANTY AND LIABILITY LIMITATIONS

PEROUSE MEDICAL warrants that this product conforms to the legal requirements, in particular those set forth in the EU Guidelines for Medical Products, 93/42/EWG. Therefore, in case this product exhibits faults that can be directly attributed to material or processing errors that occurred during its manufacture or packaging („faults“), and said faults appeared prior to the expiration date (UBD), PEROUSE MEDICAL will, at its discretion, replace the product

- a) with one of equal construction, or;
- b) with one exhibiting the identical function.

All the following conditions must be met in order to make claims within the context of this warranty:

- a) the product must not have been employed after its UBD.

- b) the fault did not initially occur after the UBD.
- c) if the product in question is a disposable product or is a product which PEROUSE MEDICAL designates as being for single use, the fault must be detected at the first use.
- d) the fault was immediately reported to PEROUSE MEDICAL or its representative.
- e) the faulty product was forwarded to PEROUSE MEDICAL or its representative immediately after the fault was reported.

This warranty in no way limits and applies in conjunction with any and all regulations governing warranties on the part of the seller and the product liability laws.

For the remainder, this warranty applies only to the extent defined above. Additional claims arising from this warranty are excluded. This warranty as well as the cited legal claims apply in place of, and to the exclusion of, all other warranties, either verbal or written, or express or implicit on the part of PEROUSE MEDICAL. None of our representatives or employees are empowered to reach differing agreements.

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In particular, PEROUSE MEDICAL assumes no liability for damages resulting from resterilization or reuse of products approved for single use only, or which PEROUSE MEDICAL designates as being for single use.

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INSTRUCTIONS FOR USE

ARTERIAL EMBOLECTOMY BALLOON CATHETER

WARNING

Read carefully instructions before use.

INDICATIONS

The arterial embolectomy balloon catheter is indicated for the removal of fresh, soft emboli and thrombi from vessels in the arterial system.

DESCRIPTION

Catheters are composed of polyamide tube with a Luer Lock hub and a tied latex balloon in its distal end. These catheters are easily identified by a color code which corresponds to their size in French. They are marked in 10 cm graduation to gauge the length of insertion and noted on the tube at the hub proximity, its naming, size and the maximal volume to inject in the balloon.

CONTRAINDICATIONS

- Do not use for patients who are allergic to latex or whose sensitivity to latex is unknown.
- Do not use in the central circulatory system
- Do not use for the removal of fibrous, adherent or calcified thrombi
- Do not use out of the arterial system
- Do not use for endarterectomy procedures
- Do not use as a vessel dilator
- Do not use in the venous system

WARNINGS AND PRECAUTIONS

- The possibility of a balloon rupture must be taken into account when considering the risks involved in any embolectomy procedures.
- A rupture can occur if the storage conditions are not respected (refer to storage part) or by contact with the atherosclerotic plaque in the vessel, particularly if it is calcified.
- To minimize the risks of vessel damage and balloon rupture, do not exceed the maximum recommended inflation volume for each size catheter (see specification table). If balloon rupture occurs, the catheter should be withdrawn.
- Do not use air to inflate the balloon: risk of air embolism in case of balloon rupture. Carbon dioxide gas is recommended for the inflation of 2F catheter (see instruction for use).
- Inspection for leaks should be performed while priming the catheter with the maximum recommended volume. Do not use if leakage is observed or inflates with difficulty.
- The medical judgment of the physician must be exercised at all time. This catheter is not made to diagnose, to watch or correct a heart or central circulation system failure.

COMPLICATIONS

As with all catheterization procedures, complications may occur. These may include local or systemic infection, local haematomas, intima disruption, arterial dissection, perforation and vessel rupture, hemorrhage, arterial thrombosis, distal embolization of blood clots and atherosclerotic plaque, air embolus, aneurysm, arterial spasm, arteriovenous fistula formation and balloon rupture with fragmentation, tip separation and distal embolization.

INSTRUCTIONS FOR USE

Catheter priming: the catheter should be primed before use.

- Select the appropriate size Luer Lock syringe for use with the catheter (see the balloon maximal volume noted next to the hub). A color code allows identifying catheters.
- Use liquid prime which is sterile and blood compatible (for example sterile physiological serum). Fill the syringe with liquid prime.
- Connect the syringe to the Luer Lock.
- Fill the catheter with the solution, remove air from the catheter by tapping it so that the air will return up into the syringe.
- Once the catheter is fully primed, the syringe is disconnected from the catheter and all residual contents are expelled from the syringe.
- The syringe is then refilled with appropriate volume of liquid and the syringe reconnected.
- The catheter is ready for use.

Note: for all inflations, use the smallest syringe capable of holding the stated maximum fluid capacity.

Inspect

- The catheter should not be used in case of leakage or if the balloon is asymmetric and inflates with difficulty.

Catheter positioning

- With the balloon deflated, insert the catheter into the vessel until and beyond the obstructing material.

Balloon inflation

- Inflate the balloon with sterile liquid or gas until it can be felt to engage the arterial wall. When no resistance is encountered, it should be assumed that the balloon has ruptured. Discontinue inflation and remove the catheter at once.

Note: the 2F catheter should be inflated with carbon dioxide gas. Gas seepage through the intact balloon will require frequent adjustment of the inflation volume.

Caution: To minimize lateral wall pressures and shear forces on the inner surface of the artery, use the smallest balloon which once inflated diameter will allow removing the obstructing material in the vessel.

Withdraw

- Remove the occlusive materials in the vessel by gently removing the catheter. During the withdrawal, it is important to adjust the balloon diameter to the varying arterial diameters by cautiously controlling the inflation volume.

PACKAGING

Sterilization method using Ethylene oxide.

Do not reuse : risk of balloon embrittlement, infection risk due to blood contamination.

Do not use the device if the packaging is damaged or opened.

Do not re-sterilize.

STERILE EO 2 LATEX 0459 2006

STORAGE

- The catheters should be stored in a dry, dark area away from heat

EXPIRY DATE

Use the catheters before expiry date indicated in the package. After expiry date, the balloon may be damaged

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August 12th, 2013

Agnes JORDAN
Quality and Regulatory
Affairs Manager

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VALIDATED

Agnes JORDAN
Corporate Regulatory
Affairs Manager

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**Additional information to the instructions for use
(maintenance documentation)
for medical devices**

**«POLYTHESE[®] IC-3GL, POLYMAILLE[®] Extra Thin,
POLYMAILLE[®] Flow and Flow Plus Heparin vascular
prostheses, balloon catheters for arterial embolectomy»**

AG

Medical devices: POLYTHESE® IC-3GL, POLYMAILLE® Extra Thin, POLYMAILLE® Flow and Flow Plus Heparin vascular prostheses, balloon catheters for arterial embolectomy

Manufacturer of POLYTHESE® IC-3GL, POLYMAILLE® Extra Thin, POLYMAILLE® Flow/Flow Plus Heparin and balloon catheters for arterial embolectomy:

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Indications:

POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses:

They are indicated to be used in cardiovascular surgery for replacement or bypass of arteries presenting aneurysm or arterial disease.

POLYMAILLE® Flow and POLYMAILLE® Flow Plus Heparin vascular prostheses:

They are indicated in cases of arterial reconstruction, primarily in the peripheral blood vessels. It can also be used in extra-anatomical reconstructions – femoro-femoral and axillofemoral.

Standardwall prostheses are also indicated for use as arterio-venous shunt prostheses in haemodialysis.

Balloon catheters for arterial embolectomy:

They are indicated for the removal of fresh, soft emboli and thrombi from vessels in the arterial system during the procedure of implantation of the vascular prostheses.

Terms of use: The above products should be only used in hospitals with special equipped facilities (the operating theater) by medical personnel with appropriate qualifications.

Technical description of the product:

POLYTHESE® IC-3GL are woven vascular prostheses, which are designed with three guide lines along the length of the prosthesis, separated from each other by 120°, and focused on the commissural line, which greatly simplifies the implantation of the prosthesis during the replacement of the aortic arch or descending part.

POLYMAILLE® Extra Thin are thin walled knitted vascular prostheses which technical characteristics completely replace a patient's artery function, up to pulse wave registration.

POLYMAILLE® FLOW and POLYMAILLE® FLOW PLUS HEPARIN are vascular prosthetic prostheses made of expanded polytetrafluoroethylene (ePTFE) of new generation. The structure of applied ePTFE allows reducing the leakage of blood at the puncture prosthesis suture. The lumen surface of the POLYMAILLE® Flow Plus Heparin prosthesis is coated with highmolecular heparin for preventing thrombosis and blood clots forming.

The following materials and substances are used for the manufacturing of POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses:

- polyester (polyethylene terephthalate - PET);
- collagen of bovine origin (1 type oxide coated), derived from the skin of Australian calves under the age of 6 months. Said material (skin) belongs to the class C (OMS classification of the European Union, according to which any possible infection is excluded). Collagen has a certificate of conformity issued by the European Pharmacopoeia number 1483, 5th edition and other processed versions of "Products with risk of BSE infection of animal origin". Certificate Number - R1 CEP 2000-343.

The following materials and substances are used for the manufacturing of POLYMAILLE® Flow and Flow Plus Heparin vascular prostheses:

2

- ePTFE;
- bioactive heparin, controlled in compliance with CPGM 7356.002 “Active pharmaceutical ingredients (API) Process inspection ” (for POLYMAILLE® Flow Plus Heparin only)

Prostheses are produced with pure undamaged surface. Vascular prostheses are flexible and durable. POLYMAILLE® Flow and Flow Plus Heparin vascular prostheses should be waterproof at a pressure of 120 mm Hg. Water permeability of POLYMAILLE® Extra Thin vascular prostheses does not exceed 5 ml/min/cm² at a pressure of 120 mm Hg. Water permeability of POLYTHESE® IC-3GL vascular prostheses does not exceed 10 ml/min/cm² at a pressure of 120 mm Hg.

Balloon catheters for arterial embolectomy are composed of polyamide tube with a Luer-Lock hub and a tied latex balloon in its distal end. Catheters are easily identified by color code which corresponds to their size in French. They are marked in 10 cm graduation to gauge the length of insertion and noted on the tube at the hub proximity, its naming, size and the maximal volume to inject in the balloon. Available sizes of catheters: from 2 to 8 Fr, length: 40, 60 and 80 cm.

Delivery set:

POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses:

- vascular prosthesis (1 piece);
- package (1 piece);
- instructions for use (1 piece);
- patient card (1 piece);
- labels for the physician (6 pieces).

POLYMAILLE® FLOW and POLYMAILLE® FLOW PLUS HEPARIN vascular prostheses:

- vascular prosthesis (1 piece);
- package (1 piece);
- instructions for use (1 piece);
- patient card (1 piece);
- labels for the physician (4 pieces).

Balloon catheters for arterial embolectomy:

- balloon catheter (1 piece);
- package (1 piece);
- instructions for use (1 piece).

Product sizes:

- for POLYTHESE® IC-3GL vascular prostheses - nominal inner diameter at rest (in mm) and effective length (in cm);
- for POLYMAILLE® Extra Thin vascular prostheses - nominal inner diameter at rest (in mm), effective length (in cm) and the length of the reinforced area (if presence);
- for POLYMAILLE® Flow and Flow Plus Heparin vascular prostheses - nominal inner diameter at rest (in mm), effective length (in cm), wall thickness (thin-walled and standard) and the presence (absence) of the reinforcement;
- for embolectomy catheters - the catheter size (in Fr) and the nominal length (in cm).

Main dimensions of the prostheses and accessories are presented in Tables 1 and 2.

Table 1.

Main dimensions of vascular prostheses

POLYTHESE® 3GL					
Catalogue number	Inner diameter (mm)	Length (cm)	Catalogue number	Inner diameter (mm)	Length (cm)
PT3GL020015	20	15	PT3GL030060	30	60
PT3GL020030	20	30	PT3GL032015	32	15
PT3GL020060	20	60	PT3GL032030	32	30
PT3GL022015	22	15	PT3GL032060	32	60
PT3GL022030	22	30	PT3GL034015	34	15
PT3GL022060	22	60	PT3GL034030	34	30
PT3GL024015	24	15	PT3GL034060	34	60
PT3GL024030	24	30	PT3GL036015	36	15
PT3GL024060	24	60	PT3GL036030	36	30
PT3GL026015	26	15	PT3GL036060	36	60
PT3GL026030	26	30	PT3GL038015	38	15
PT3GL026060	26	60	PT3GL038030	38	30
PT3GL028015	28	15	PT3GL038060	38	60
PT3GL028030	28	30			

POLYMAILLE® Extra Thin							
Catalogue number	Inner diameter (mm)	Length (cm)	Length of the reinforced area (cm)	Catalogue number	Inner diameter (mm)	Length (cm)	Length of the reinforced area (cm)
PMET060040	6	40	-	PMETR06100060	6	100	60
PMET060070	6	70	-	PMETR07040015	7	40	15
PMET060100	6	100	-	PMETR07040020	7	40	20
PMET070040	7	40	-	PMETR07070020	7	70	20
PMET070070	7	70	-	PMETR07070040	7	70	40
PMET070100	7	100	-	PMETR07100060	7	100	60
PMET080040	8	40	-	PMETR08040015	8	40	15
PMET080070	8	70	-	PMETR08040020	8	40	20
PMET080100	8	100	-	PMETR08070020	8	70	20
PMET100040	10	40	-	PMETR08070040	8	70	40
PMET100070	10	70	-	PMETR08100060	8	100	60
PMET100100	10	100	-	PMETR10040015	10	40	15
PMETR06040015	6	40	15	PMETR10040020	10	40	20
PMETR06040020	6	40	20	PMETR10070020	10	70	20
PMETR06070020	6	70	20	PMETR10070040	10	70	40
PMETR06070040	6	70	40	PMETR10100060	10	100	60

A

POLYMAILLE® Flow		
Inner diameter (mm)	Length (cm)	Catalogue number Thin-walled
6	20	PFT0620
7	20	PFT0720
8	20	PFT0820
6	40	PFT0640
7	40	PFT0740
8	40	PFT0840
6	50	PFT0650
7	50	PFT0750
8	50	PFT0850
6	80	PFT0680
7	80	PFT0780
8	80	PFT0880

POLYMAILLE® Flow		
Inner diameter (mm)	Length (cm)	Catalogue number Thin-walled, reinforced
6	40	PFTR0640
7	40	PFTR0740
8	40	PFTR0840
6	50	PFTR0650
7	50	PFTR0750
8	50	PFTR0850
6	80	PFTR0680
7	80	PFTR0780
8	80	PFTR0880

POLYMAILLE® Flow		
Inner diameter (mm)	Length (cm)	Catalogue number Standard wall
6	10	PFS0610HD
7	10	PFS0710HD
6	20	PFS0620HD
7	20	PFS0720HD
6	40	PFS0640
7	40	PFS0740
8	40	PFS0840
6	50	PFS0650
7	50	PFS0750
8	50	PFS0850
6	80	PFS0680
7	80	PFS0780
8	80	PFS0880

5

POLYMAILLE® Flow		
Inner diameter (mm)	Length (cm)	Catalogue number Standard wall, reinforced
6	40	PFSR0640
7	40	PFSR0740
8	40	PFSR0840
6	50	PFSR0650
7	50	PFSR0750
8	50	PFSR0850
6	80	PFSR0680
7	80	PFSR0780
8	80	PFSR0880

POLYMAILLE® Flow Plus Heparin		
Inner diameter (mm)	Length (cm)	Catalogue number Thin-walled
6	20	PFHT0620
7	20	PFHT0720
8	20	PFHT0820
5	40	PFHT0540
6	40	PFHT0640
7	40	PFHT0740
8	40	PFHT0840
5	50	PFHT0550
6	50	PFHT0650
7	50	PFHT0750
8	50	PFHT0850
5	80	PFHT0580
6	80	PFHT0680
7	80	PFHT0780
8	80	PFHT0880

POLYMAILLE® Flow Plus Heparin		
Inner diameter (mm)	Length (cm)	Catalogue number Standard wall
6	10	PFHS0610HD
7	10	PFHS0710HD
6	20	PFHS0620HD
7	20	PFHS0720HD
5	40	PFHS0540
6	40	PFHS0640
7	40	PFHS0740
8	40	PFHS0840
5	50	PFHS0550
6	50	PFHS0650
7	50	PFHS0750
8	50	PFHS0850

Table 2.

Main dimensions of balloon catheters for arterial embolectomy

Balloon catheters for arterial embolectomy					
Catalogue number	Size (Fr)	Length (cm)	Color	Maximum fluid capacity	Inflated balloon diameter (mm)
CEM802	2	80		0.2 (gas)	4
CEM602	2	60		0.2 (gas)	4
CEM402	2	40		0.2 (gas)	4
CEM803	3	80		0.20 ml	5
CEM603	3	60		0.20 ml	5
CEM403	3	40		0.20 ml	5
CEM804	4	80		0.75 ml	7
CEM604	4	60		0.75 ml	7
CEM404	4	40		0.75 ml	7
CEM805	5	80		1.50 ml	9
CEM605	5	60		1.50 ml	9
CEM405	5	40		1.50 ml	9
CEM806	6	80		2.00 ml	11
CEM606	6	60		2.00 ml	11
CEM406	6	40		2.00 ml	11
CEM807	7	80		2.50 ml	13
CEM407	7	40		2.50 ml	13
CEM808	8	80		3.00 ml	14
CEM608	8	60		3.00 ml	14

Sterilization conditions: Vascular prostheses, balloon catheters for arterial embolectomy are supplied sterile.

Sterilization method of POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses is β -radiation.

POLYMAILLE® Flow and Flow Plus Heparin vascular prostheses, balloon catheters for arterial embolectomy are sterilized using Ethylene oxide.

Disposable medical products. Do not re-sterilized.

Vascular prostheses and balloon catheters for arterial embolectomy are non-cytotoxic.

Each serie of products produced by PEROUSE MEDICAL are inspected to be in compliance with factory standards, sterilization process validation and receives the certificate of compliance/non-compliance with existing standards. Products, which are not compliant with the standards requirements, cannot be further distributed.

Packaging: Each product is individually packaged.

Vascular prostheses are placed in a double blister and then packaged in a box. POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses are also placed in an aluminum pouch prior to packaging in box.

Balloon catheters for arterial embolectomy are packaged in a rigid cane then in a peel-away pouch. One side of the pouch is made of transparent material which gives the opportunity to inspect the product during its preparation for surgery.

Documentation supplied with the product is also placed in the box.

Transportation and storage:

POLYTHESE® IC-3GL and POLYMAILLE® Extra Thin vascular prostheses should be transported and stored away from light, cold, heat at temperatures between 10°C and 35°C and relative humidity of 80%.

POLYMAILLE® Flow should be transported and stored in a cool, clean and dry place.

POLYMAILLE® Flow Plus Heparin should be transported and stored in a cool, clean, dry environment at temperatures between 10°C and 40°C.

Balloon catheters for arterial embolectomy should be transported and stored in a dry, dark area away from heat.

Shelf life of POLYTHESE® IC-3GL, POLYMAILLE® Extra Thin, POLYMAILLE® Flow vascular prostheses is 5 years, POLYMAILLE® Flow Plus Heparin vascular prostheses is 2 years, balloon catheters for arterial embolectomy is 2.5 years.

Recommendations:

Vascular prostheses, balloon catheters for arterial embolectomy should not be used after the expiration date and/or if packaging is damaged.

Balloon catheters for arterial embolectomy should not be used in case of leakage or if the balloon is asymmetric and inflates with difficulty.

In all the above situations, and in case the inability to use the product according to its intended purpose, medical device must be disposed immediately.

Use of POLYTHESE® IC-3GL, POLYMAILLE® Extra Thin, Flow and Flow Plus Heparin vascular prostheses, balloon catheters for arterial embolectomy must be made in accordance with applicable health rules and regulations of the Russian Federation (SanPiN) for the handling of medical waste and the Guidelines for disinfection, and sterilization of medical devices (MU).

Warranty and liability limitations:

PEROUSE MEDICAL warrants that these products are in compliance with the legal requirements, in particular those set forth in the EU Directive 93/42/EEC concerning medical devices. Therefore, in case the product exhibits non-conformities or defects that can be directly attributed to material or processing errors that occurred during its manufacture or packaging ("non-conformities or defects"), and said non-conformities or defects appeared prior to the expiration date (UBD), PEROUSE MEDICAL will, at its discretion, replace the product:

- with one of equal construction, or
- with one exhibiting the identical function.

All the following conditions must have been met in order to make claims within the context of this warranty:

- the product must not have been used after its UBD.
- The non-conformity or defect did not initially occur after the UBD.
- If the product in question is a disposable product or is a product, which PEROUSE MEDICAL designates as being for single use, the fault must be detected at the first use.
- the non-conformity or defect was immediately reported to PEROUSE MEDICAL or its representative.
- the non-compliant or defective product was forwarded to PEROUSE MEDICAL or its representative immediately after the non-conformity or defect was reported.

The above warranty in no way limits and applies in conjunction with any and all regulations governing warranties on the part of the seller and the product liability laws.

For the remainder, this warranty applies only to the extent defined above. Additional claims arising from this warranty are excluded. This warranty, as well as the cited legal claims applies in place of, and to the exclusion of, all other warranties, either verbal or written, or express or implicit in the part of PEROUSE MEDICAL. None of our representatives or employees is empowered to reach differing agreements.

Due to numerous factors such as shipping, storage, and handling on the part of the user, all of which lie outside the control of PEROUSE MEDICAL, PEROUSE MEDICAL is not liable, in particular, for damages of any type resulting directly or indirectly from the product or its employment unless said damages can be shown to be the result of a failure on the part of the product.

In particular, PEROUSE MEDICAL assumes no liability for damages resulting from sterilization or reuse of products approved for single use, or which PEROUSE MEDICAL designates as being for single use.

Instructions for use must be accurately studied before using of the product.

Only highly qualified medical staff with special training to carry out the relevant procedures should operate medical devices manufactured by PEROUSE MEDICAL.

All the pending disputes should be resolved in the courts at the location of PEROUSE MEDICAL, in accordance with French law.

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
1. Expiry date:

POLYTHESE® IC 3GL and POLYMAILLE® EXTRA THIN vascular prostheses are sterilized with β radiation. Expiry date for these products is 5 years.

2. Storage conditions:

The products should be stored away from light, cold, heat ($10^{\circ}\text{C} < T < 35^{\circ}\text{C}$) and humidity ($\text{RH} < 80\%$).

3. Sterilization

The validation process was conducted by IONISOS at Chaumesnil (France), according to NF EN ISO 11137-1:2006 standard (sterilization of healthcare -Radiation - Part 1: Requirements for development, validation and routine control of a sterilization process for medical devices). The scope of the validation process is to confirm the ability of the radiation sterilization cycle to deliver and maintain a sterility assurance level (SAL) of 10^{-6} .

3.1 Sterilization conditions:

Minimum radiation dose: 25 kGy

Maximum radiation dose: 50 kGy

3.2 Sterilization routine controls:

PERFORMED TESTS	SPECIFICATION	FREQUENCY	APPLIED STANDARD
Control of the sterilization cycle (β radiation)	In compliance with the validated cycle	At each sterilization	NF EN ISO 11137-1
Estimation of the population of micro-organisms on product	Acceptable < 300 CFU*/article	1 non sterile device by product family at each sterilization	NF EN ISO 11737-1
Pyrogen test	European Pharmacopeia § 2.6.8	1 test once a month by group family	European Pharmacopeia

*CFU: Colony Forming Unit

3.3 Sterilization effectiveness control

Sterilization dose audits (by alternating family product references) are regularly performed according to the frequencies defined in paragraph 12.1.3.1 of NF EN ISO 11137-1 standard (2006).

Ivry Le Temple,
March 6th 2014

Marie Noelle EROUT
Quality Manager-Ivry site

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POLYMAILLE FLOW – Storage and Sterilization

1. Expiry date

Products are sterilized with ethylene oxide. Expiry date for POLYMAILLE FLOW vascular prosthesis is 5 years.

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Aeration Room	Area [m ²]	High [m]	Room Volume [m ³]	Rate of Air Extraction [1/hour]	Volume of extracted air [m ³ /hour]	Room changes per hour	Temperature [°C]
Aeration Room I	428.6	10.5	4.500	1.5	6.750	10	35 - 45
Manufacturer: InfraServ Wiesbaden, Year of installation 2005							

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Parameter	Unit	Set-point and specifications
		Cycle #15 Rev. 1
Preconditioning		
Temperature	°C	40 ± 5
RH	%	60 ± 15
Time	hh:mm	16:00 – 72:00, ≤ 16:00
Transfer time into the chamber	min	≤ 60; ≥ 60
Chamber Cycle (recirculation blower in use)		
Evacuation Pressure	mbarA	70 (± 5)
Leak test	hh:mm	00:10
Leak test maximum leak rate	mbar/min	1.5
Nitrogen dilution (sub phase: steam injection)	mbar	20 (± 5)
Nitrogen dilution (sub phase: nitrogen injection)	mbarA	250 (± 5)
Nitrogen dilution (sub phase: re-evacuation)	mbarA	70 (± 4)
Nitrogen dilution number of pulses	-	2
Steam injection 1 (Δ p)	mbar	30 (± 5)
Steam dwell (static)	hh:mm	00:30
Steam dwell end pressure	mbarA	100 (± 5)
ETO injection 1	mbarA	430 (± 5)
Nitrogen top up 1	mbarA	455 (± 5)
ETO injection 2	mbarA	555 (± 5); 545 (± 5)
Nitrogen top up 2	mbarA	610 (± 5)
Gas dwell pressure range	mbarA	593 - 625
ETO concentration (theoretical, based on pressure differential)	mg/l	717
ETO dwell time	hh:mm	05:00, 02:30
ETO dwell temperature Set point (average of 3 Pt-100s)	°C	45; 43
Max Range Chamber Temp during dwell time	°C	40 – 50; 38 - 48
Evacuation pressure	mbarA	80 (± 10)
Number of wash cycles (nitrogen)	-	2
Nitrogen washes: injection pressure	mbarA	880 (± 10)
Nitrogen washes: evacuation pressure	mbarA	90 (± 10)
Number of wash cycles (air)	-	2
Air washes: injection pressure	mbarA	880 (± 10)
Air washes: evacuation pressure	mbarA	100 (± 10)
Final release	Atmospheric pressure	
Aeration		
Temperature	°C	40 ± 5

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POLYMAILLE FLOW PLUS HEPARIN – Storage and Sterilization

1. Expiry date

Products are sterilized with ethylene oxide. Expiry date for POLYMAILLE FLOW PLUS HEPARIN vascular prosthesis is 2 years.

2. Storage conditions

The products should be stored in a dry and cool place.

3. Sterilization

The validation process was conducted by Sterigenics, Wiesbaden, according to ISO 11135-1:2007 standard by using the overkill method in 2 half cycles and 1 full cycle. The scope of the validation process is to confirm the ability of the EO sterilization cycle to deliver and maintain a sterility assurance level (SAL) of 10^{-6} .

3.1. Specifications

Preconditioning

Specification of Sterigenics Wiesbaden Preconditioning:

Preconditioning (PC)	Area [m²]	High [m]	Room Volume [m³]	Rate of Air Extraction [1/hour]	Room changes per hour	Relative Humidity [%]	Temperature [°C]
PC B	44.0	3.8	167	0	20	45 - 75	35 – 45
PC D	174.0	10.0	1740				
Manufacturer: InfraServ Wiesbaden, Year of installation 2005							

Sterilization (exposure to EO)

Specification of Sterigenics Wiesbaden Chamber B:

Chamber code	Date of Manufacture	Number of pallets max.	Maximum height of pallets	Usable volume (usv - m ³)	Chamber heat
B	1994	13 Euro	2.0 m	24.96 m ³	Floor, ceiling, door, walls
Manufacturer: Sterichem France/Getinge Sweden					

Aeration

Specification of Sterigenics Wiesbaden Aeration Room:

Aeration Room	Area [m ²]	High [m]	Room Volume [m ³]	Rate of Air Extraction [1/hour]	Volume of extracted air [m ³ /hour]	Room changes per hour	Temperature [°C]
Aeration Room I	428.6	10.5	4.500	1.5	6.750	10	35 – 45
Manufacturer: InfraServ Wiesbaden, Year of installation 2005							

3.2. Ethylene oxide

All validation runs were conducted using 99.9 % ethylene oxide, which is certified for each run separately in Annexes of each cycle in the validation report. Ethylene oxide used for this study was delivered by Sterigenics' approved supplier: Gerling Holz (Hanau, Germany).

3.3. Cycle parameters

The following table showed the cycle parameters for a full cycle and its derivatives. For microbiological performance qualification, certain parameters were reduced relative to the proposed routine process, marked in red for the half cycle.

Parameter	Unit	Set-point and specifications
		Cycle #15 Rev. 1
Preconditioning		
Temperature	°C	40 ± 5
RH	%	60 ± 15
Time	hh:mm	16:00 – 72:00, ≤ 16:00
Transfer time into the chamber	min	≤ 60; ≥ 60
Chamber Cycle (recirculation blower in use)		
Evacuation Pressure	mbarA	70 (± 5)
Leak test	hh:mm	00:10
Leak test maximum leak rate	mbar/min	1.5
Nitrogen dilution (sub phase: steam injection)	mbar	20 (± 5)
Nitrogen dilution (sub phase: nitrogen injection)	mbarA	250 (± 5)
Nitrogen dilution (sub phase: re-evacuation)	mbarA	70 (± 4)
Nitrogen dilution number of pulses	-	2
Steam injection 1 (Δ p)	mbar	30 (± 5)
Steam dwell (static)	hh:mm	00:30
Steam dwell end pressure	mbarA	100 (± 5)
ETO injection 1	mbarA	430 (± 5)
Nitrogen top up 1	mbarA	455 (± 5)
ETO injection 2	mbarA	555 (± 5); 545 (± 5)
Nitrogen top up 2	mbarA	610 (± 5)
Gas dwell pressure range	mbarA	593 - 625
ETO concentration (theoretical, based on pressure differential)	mg/l	717
ETO dwell time	hh:mm	05:00, 02:30
ETO dwell temperature Set point (average of 3 Pt-100s)	°C	45; 48
Max Range Chamber Temp during dwell time	°C	40 – 50; 38 - 48
Evacuation pressure	mbarA	80 (± 10)
Number of wash cycles (nitrogen)	-	2
Nitrogen washes: injection pressure	mbarA	880 (± 10)
Nitrogen washes: evacuation pressure	mbarA	90 (± 10)
Number of wash cycles (air)	-	2
Air washes: injection pressure	mbarA	880 (± 10)
Air washes: evacuation pressure	mbarA	100 (± 10)
Final release	Atmospheric pressure	
Aeration		
Temperature	°C	40 ± 5

3.4. Residual Testing

An analysis on residual Ethylene Oxide (ETO) and Ethylene Chlorohydrine (ECH) was performed on a worst case product.


The samples for residual testing were processed with one full cycle. They were placed in the aeration room. 1 sample of each product type was removed from the load after 0, 48 and 96 hours of aeration. The products were frozen at -20 °C until the analysis was performed.

A minimum of 4 days of heated aeration is required to fulfil the requirements of ISO 10993-7:2008.

4. Routine controls

Bioburden / LAL testing: Products are monitored at regular intervals in routine controls.

Ivry Le Temple,
March 6th 2014

Marie-Noelle EROUT
Quality Manager – Ivry Site

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Balloon catheters for arterial embolectomy- Storage and conservation

EXPIRY DATE:

Balloon catheters for arterial embolectomy are sterilized with Ethylene Oxide. Expiry date for these products is two and half years.

STORAGE CONDITIONS:

The products should be stored away from light and heat.

STERILIZATION

It is carried out in accordance with NF EN ISO 11135-1 standard.

Sterilization conditions:

PRECONDITIONNING:

Temperature before preconditioning: 17°C minimum

Chamber n°: C, E

Registered temperature: 44-53°C

Registered humidity: 48-70% RH

Duration: 11h30-72h00

TRANSFERT DURATION: 45 min maximum

STERILIZATION:

Chamber n°: C

Check value: -34801

Initial vacuum: 90-110 mbar A

Inert number: 1 minimum

ΔP vapour injection: 52-60 mbar

Nitrogen injection pressure: 590-610 mbar A

Evacuation pressure: 90-110 mbar A

ΔP vapour injection: 52-60 mbar A

Vapour exposure duration : 20-30 min

Vapour exposure pressure : 144-165 mbar A

N₂ first injection : 185-205 mbar A

OE injection end : 670-690 mbar A

ΔP OE injection: 465-505 mbar

N₂ second injection: 790-810 mbar A

OE mass: 15 kg minimum

Gas exposure duration: 02h59-04h00

Gas exposure temperature: 45-53°C

Gas exposure pressure: 619-900 mbar A

Evacuation pressure: 90-110 mbar A

Rinse number: 3 minimum

Nitrogen injection pressure: 490-510 mbar A

Evacuation pressure: 90-110 mbar A

Cycle end pressure: 890 mbar A minimum

The sterilization process was described in the validation report N°RVSR271110 rev00 and the kinetic desorption in the validation report N°RVAR150312 rev00.

The routine sterilization is released according to:

- The sterilization process parameters,
- The internal and external biological indicators,
- The certificate of desorption of OE,
- The endotoxin test for each batch.

Ivry Le Temple,
March 21st, 2014

Marie Noelle EROUT
Plant Quality Manager

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