GE Healthcare

T2100-ST1

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

Notary Certificate

State of Kansas
County of Harvey

Before me, the undersigned notary public, this day, personally, appeared to me known or proven, who being duly sworn according to the law, deposes the following:

I, Krishyn Caldwell, hereby certify that the attached document is
☒ an original document
☐ a true and correct copy of an original document


to the best of my knowledge and belief.

Krishyn L Caldwell
Signature of notarial officer, Notary Public

04-19-2024
Date

My commission expires: 08-22-2026

Notary Stamp

Publication Information

The information in this manual applies only to the T2100-ST1 (GE PN 2097357-001) and T2100-ST2 (GE PN 2097357-002) treadmills. It does not apply to earlier versions. Due to continuing product innovation, specifications in this manual are subject to change without notice.

T2100-ST1, T2100-ST2, CASE, CardioSoft, and MAC are trademarks owned by GE Medical Systems Information Technologies, Inc., a General Electric Company going to market as GE Healthcare. All other marks are the properties of their respective owners.

This product complies with the regulatory requirements concerning medical devices from the following bodies:

The T2100-ST1 and T2100-ST2 treadmills meet the following safety and regulatory standards for FDA Class 1 motor operated physical medicine machines. They have been tested by Intertek Testing Services N.A Inc. and are listed by Engineering Testing Laboratories (ETL). However, the ultimate conformance to IEC 6060-1 2005-3rd edition is the responsibility of the system integrator when combined with other equipment. Additionally, all motorized equipment is potentially dangerous if used incorrectly. Before using the T2100-ST1 or T2100-ST2 treadmill, follow all precautions listed in this manual and read the entire operator's manual thoroughly. Use the T2100-ST1 and T2100-ST2 treadmills only as described.

Revision History

The document part number and revision appear at the bottom of each page. The revision identifies the document's update level. The revision history of this document is summarized in the following table.

Revision	Publication Date	Description
A	20 August 2016	Internal release
B	9 September 2016	Initial public release
C	18 November 2016	Power cord configuration
D	4 April 2017	Updated eStop button
E	12 June 2017	Added power cord configuration
F	27 June 2018	Revised regulatory & safety conformance
G	22 January 2019	Added EMC Declaration and 4 th Edition Configuration
H	14 December 2019	Update Speed & Calibration Procedure
J	02 March 2021	Added MDR Labeling, Uncrating Instructions & Reprocessing Instructions
K	05 December 2021	Implementation of New Motor and Drive Combination, 3 Wire Hall Effect Speed Sensor and Harness, Magnetic Tether Switch with STOP Behavior Information, Calibration Procedure Software Update and Firmware Change.
L	22 June 2022	update of 11/32" Wrench and Socket to required tools and added a Note on the 11-wire cable.
M	13 March 2023	Implemented alternative power switch and T2100-ST Tether Guard

Revision	Publication Date	Description
N	7 July 2023	Updated the EC REP address, added mounting and dismounting method steps, added additional warnings in regard to the Emergency Stop Button and Patient Safety Tether functionality and behavior, removed contradicting warnings, updated text for the Patient Safety Tether, removed reference to T2100®, updated 16 hour to minimum of 4 hours functional test. Updated replaced part numbers. Updated. Updated to latest firmware on the Smart Power Supply Board. Updated latest Calibration software.

To access other GE Healthcare Diagnostic Cardiology documents, go to the Common Documentation Library (CDL), located at www.gehealthcare.com/documents, and click Cardiology.

Contents

1	Introduction	11
	Intended User	11
	Indications for Use	11
	Prescription Device Statement.....	11
	Regulatory and Safety Information	11
	Safety Conventions	12
	Safety Hazards	12
	Classification of Medical Device	15
	Regulatory and Safety Conformance	15
	Responsibility of the Manufacturer	19
	Responsibility of the Customer	20
	Product and Package Information	22
	Symbols	22
	Equipment Identification.....	30
	Product Label.....	30
	Service Information.....	31
	Service Requirements	31
	Warranty Information	31
	Additional Assistance	31
	Reporting of serious incidents.....	32
	Manual Information	32
	Manual Purpose.....	32
	Document Conventions.....	32
	Related Documents.....	34
	Training	36
2	Product Overview	37
	Safety Systems	38
	Treadmill.....	38
	Drive System.....	38
	Speed Range	38
	Incline Range.....	38
	Running Surface	38
	Communication Ports	39
	Floor Surface Footprint.....	39
	Operating and Storage Condition Recommendations.....	39
	Power Requirements	39

3	Assembly and Setup.....	43
	Safe Handling Guidelines.....	43
	Uncrating Instructions.....	44
	Initial Setup.....	47
	Location.....	48
	Electrical Safety Tests.....	49
	AC Line Voltage Test.....	49
	Power Cord and Plug.....	50
	Ground Continuity Test.....	51
	Conducting Leakage Tests.....	51
	Leakage Test Diagrams.....	52
	Secure Cables.....	53
	Running Belt Tracking Adjustment.....	53
	Running Belt Tension Adjustment.....	53
	Drive Belt Tension Adjustment.....	54
	Test Plug Procedure.....	54
	Communication Ports.....	55
	Using the test plug.....	55
	Setup and Connection to Host.....	57
	CASE connection and setup:.....	57
	CardioSoft/CS connection and setup:.....	59
	MAC 5500 ST connection and setup:.....	65
	MAC 2000 ST connection and setup:.....	66
4	Preventive Maintenance.....	69
	Reprocessing Instructions.....	69
	Daily Maintenance.....	70
	Weekly Maintenance.....	70
	Monthly Maintenance.....	70
	Semiannual Maintenance.....	70
	Belt Cleaning and Inspection.....	71
	Running Belt Tracking Adjustment.....	71
	Running Belt Tension Adjustment.....	72
	Drive Belt Tension Adjustment.....	73
	Exterior Care.....	74
	Elevation Screw Lubrication.....	75
	Running Deck Maintenance.....	75
5	Theory of Operation.....	76

Scope	76
T2100-ST Series Block Diagram.....	77
Smart Power Supply (PCB) Overview	79
Software Requirements.....	81
Speed Control	81
Elevation Control	82
ESTOP.....	82
Communication Hardware (RS-232 Configuration).....	83
Communication Hardware (USB Configuration)	83
Self-Test Mode.....	84
Calibration	84
DIP Switch Settings.....	85
Electrical Inputs.....	85
Electrical Outputs.....	86
Electrical Connections.....	86
Physical Requirements and Restrictions	87
6 Troubleshooting	88
Troubleshooting Guidance Base on Error Code	89
Incoming Power 110-240VAC Flow Chart 1A	91
Incoming Power Inline Filter Flow Chart 1B	92
Incoming Power Drive PC2303-012-N Flow Chart 1C	93
Smart Power Supply Incoming Power 110-240VAC Flow Chart 1D	94
Stop Tether Switch Flow Chart 1F.....	96
Communication RS232 Flow Chart 1G	97
Communication USB Flow Chart 1H.....	98
Smart Power Supply Error Code Identification Flow Chart 1I.....	99
Smart Power Supply Error Code 1 Flow Chart 1J "Bad Calibration Error (1)" 100	
Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)" 101	
Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)"	102
Smart Power Supply Error Code 3 Flow Chart 1L "Zero Speed Motor Controller Error (3)"	103
Smart Power Supply Error Code 4 Flow Chart 1M "Over Speed Motor Controller Error (4)"	104

Smart Power Supply Error Code 5 Flow Chart 1N "Zero Speed External Sensor Error (5)"	105
Smart Power Supply Error Code 6 Flow Chart 1O "Over Speed External Sensor Error (6)"	106
Smart Power Supply Error Code 7 Flow Chart 1P "Speed Compare Error (7)"	107
Smart Power Supply Error Code 8 Flow Chart 1L "Motor Controller Fault Signal (8)"	109
Smart Power Supply Error Code 9 Flow Chart 1R "Belt Start Reject Error (9)"	110
Kollmorgen Drive PC2303-012-N status LED CODE list	111
Quantum Drive status CODE list.....	111
Running Belt High Speed Application	112
System Flash Log Retrieval.....	112
7 Field Replaceable Units.....	115
Final Assembly.....	116
Final Assembly Circuit Board Connection	116
Deck Assembly	120
Frame Assembly	122
Motor Pan Assembly Kollmorgen.....	124
Motor Pan Assembly Kollmorgen Wiring	127
Motor Pan Assembly Quantum/TruTech	131
Motor Pan Assembly Quantum/TruTech Wiring	134
Motor Assembly Kollmorgen.....	138
Motor Assembly Quantum/TruTech.....	138
Drive Board Assembly Quantum/TruTech	139
Circuit Board Assembly FGLF0495-1.....	141
Circuit Board Assembly FGLF0495-3.....	142
Stop Tether Switch Assembly	144
Emergency Stop Assembly	145
Communication Cables.....	145
Calibration Software	146
Firmware Upgrade	146
Optional Accessories	146
8 Component/FRU Removal and Replacement	147
Established Component Replacement Times.....	147

Comprehensive Tool List.....	148
Replacing the Hood.....	149
Replacing the Kollmorgen Drive Motor	150
Replacing the Quantum/TruTech Drive Motor Assy	153
Replacing the Kollmorgen Drive Board Assy	155
Replacing the Quantum/TruTech Drive Board Assy	157
Replacing Kollmorgen Motor & Drive Board with Quantum/TruTech Motor & Drive Board Using Conversion Kit.....	159
Replacing the Front Roller	164
Replacing the Rear Roller	165
Replacing the Running Belt	167
Replacing the Running Deck.....	170
Replacing the Deck Cushion.....	171
Replacing the Motor Drive Belt	173
Replacing/Adjusting the Elevation Actuator	175
Replacing the Speed Sensor from 2 Wire to 3 Wire	176
Replacing the Speed Sensor 3 Wire to 3 Wire.....	178
Replacing the Smart Power Supply Relay Board	179
Replacing the Circuit Breaker	181
Relocating the Emergency Stop Button (ESB) Assy and Stop Tether Switch (STS) Assy.....	182
Replacing the Emergency Stop Button (ESB) Assy	184
Replacing the Stop Tether Switch (STS) Assy	186
Replacing the Right or Left Handrail	187
Replacing the Center Handrail.....	189
Removing and Reinstalling the Handrails for Moving	190
Replacing the Main Power Switch	191
Adjusting the Run Belt Tracking and Tension.....	194
Adjusting belt tracking.....	194
Adjusting belt tension.....	195
Replacing the Power cord	195
Replacing the CE filter	196
Replacing the Relay.....	197
9 Calibration.....	199
Installing the T2100-ST Series Calibration Software	199
Installing the Serial Drivers	201

Software Interface Instruction T2100-ST Series Calibration.....	205
Calibrating Speed.....	206
Calibrating Elevation.....	209
Elevation Chart.....	212
Elevation Actuator	213
Calibration Software Tabs.....	213
Basic Control Tab	213
Advanced Control Tab	214
Speed Cal Tab	216
Grade Cal Tab	217
Error Code Tab.....	217
System Log Tab.....	218
Quantum Tab	218
10 Functional / Post Repair / Preventative Maintenance Checkout Procedures	221
Tools Required	221
Checkout Procedures.....	221
Visual/Functional Checks	222
Operational Checks.....	222
Speed Calibration.....	222
Elevation Calibration.....	223
Test Plug Check.....	223
Walking Belt Tension / Tracking	223
Walking Belt Speed Verification.....	223
Host Communication.....	224
Electrical Safety Checks	224
Preventative Maintenance Visual/Functional Checks	225
Procedural Lists.....	225
FRU Repairs / Exchange.....	226
Preventative Maintenance Visual/Functional Checks	228
A Maintenance Log	229
Reference Documents	231
11 Reprocessing Instructions	232

1

Introduction

This document describes the T2100-ST1 and T2100-ST2 treadmills also referred to as the “system”, “device”, or “product”. The document is intended to be used by service personnel.

This chapter provides general information required for the proper use of the system and this manual. Familiarize yourself with this information before using the system.

Intended User

This manual is intended for service personnel responsible for the maintenance and repair of the T2100-ST1 and T2100-ST2 treadmills.

Indications for Use

The T2100-ST1 and T2100-ST2 treadmills are designed for cardiac stress testing.

Prescription Device Statement

CAUTION:

United States federal law restricts this device to sale by, or on the order of, a physician.

Regulatory and Safety Information

This section provides information about the safe use and regulatory compliance of this system. Familiarize yourself with this information, read and understand all instructions before attempting to use this system. The system was designed and manufactured to the appropriate medical regulations and controls.

NOTE:

Disregarding the safety information provided in this manual is considered abnormal use of this system and could result in injury, loss of data, and void any existing product warranties.

Introduction

Safety Conventions

A **Hazard** is a source of potential injury to a person, property, or the system.

This manual uses the terms DANGER, WARNING, CAUTION, and NOTICE to point out hazards and to designate a degree or level of seriousness. Familiarize yourself with the following definitions and their significance.

Definition of Safety Conventions

Convention	Definition
DANGER	Indicates an imminent hazard, which, if not avoided, will result in death or serious injury.
WARNING	Indicates a potential hazard or unsafe practice, which, if not avoided, could result in death or serious injury.
CAUTION	Indicates a potential hazard or unsafe practice, which, if not avoided, could result in moderate or minor injury.
NOTICE	Indicates a potential hazard or unsafe practice, which, if not avoided, could result in loss or destruction of property or data.

Safety Hazards

The following messages are applicable for the system. Specific messages may also appear elsewhere in the manual.

WARNING:

The T2100-ST1 and T2100-ST2 treadmills are manufactured to exacting standards both in physical form and in component selection. The components used in our products have been selected with performance and medical safety in mind. The treadmills have been engineered and certified to conform to the list of medical and safety regulatory standards which appear on the next page. Modification or part substitution of any kind is strictly forbidden. Any deviation in component replacement, physical or electrical modification will result in loss of medical safety certification and warranty of this product. Modifications to this equipment may put the patient at risk of electrical shock or hardware malfunction.

Contact GE Healthcare Service department for all your repair part needs.

WARNING:

The T2100-ST1 and T2100-ST2 treadmills must be grounded to reduce the risk of electrical shock. If a malfunction occurs, earth grounding provides a path of least resistance for an electric current. Ungrounded connections must not be used.

No other equipment may be used on the electrical circuit with the treadmills. Do not use extension cords. Using a shared or unreliable circuit can also cause the treadmills to

Introduction

unexpectedly shut off, potentially resulting in injury to the patient.

Ensure the master power switch is in the off position before plugging in the T2100-ST1 or T2100-ST2. A power surge could damage the sophisticated electronic system of the treadmills.

WARNING:

Before permitting anyone to use the T2100-ST Series, do the following:

- Warn each user about the risk of falling while the belt is in motion.
- Stress the need for caution.
- Wait until the treadmill belt is moving before stepping onto the belt.
- Demonstrate the proper mounting and dismounting methods.
- Mounting methods:

Step 1: Instruct patient to straddle the belt and hold on to the handrails during initialization of the treadmill before starting the test

Step 2: Send start/run command to the treadmill with patient still straddling the belt, wait 30 seconds to ensure there is no uncommanded motion

Step 3: Instruct patient to carefully step on the belt while holding on to the handrails, begin walking and continue with the test

- Dismounting methods:

Step 1: Slow the running belt to its minimum speed, then stop the belt

Step 2: Instruct patient to carefully step off the belt while holding on to the handrails

-
- Show each user how to use the T2100-ST Series as described in this manual.
- Ask each user to perform a supervised "test usage" at minimum belt speed to review and practice usage techniques.
- Observe all the precautions listed under "Customer "on page 20 to reduce the possibility of serious injury as a result of falling or a loss of balance.

WARNING:

Introduction

Serious injury or death could result from electrical shock. To reduce the possibility of electrical shock, carefully observe the following precautions.

- To disconnect the treadmill, set the power switch to the OFF position, and remove the plug from the outlet. When the power is off, the green light on the power switch is dark.
- Never operate the unit with a damaged power cord or plug.
- Power cord should be routed through frame-mounted clamp and kept clear of the elevation mechanism.
- Communication cable should be routed through the communication cable clamp and kept clear of the elevation mechanism.
- Keep the power cord and communication cable out of traffic areas and away from heated surfaces.
- Never use extension cords.
- Never operate the unit when it is wet.
- Never operate the unit if it is not operating properly.
- Always unplug the machine before service or maintenance is performed.
- Treadmill should be serviced by authorized technicians only.
- Operator should report any electrical shock when touching the treadmill and discontinue use immediately.
- Never use the treadmill outdoors.
- Immediately discontinue use and unplug the treadmill if you smell the distinctive odor of hot electrical components.

WARNING:

Serious injury or death could result from electrical shock occurring during defibrillation. Never allow patient or operators near treadmill during defibrillation.

WARNING:

Consult your physician prior to using this appliance to determine the patient's physical readiness and capabilities. Stop exercising immediately and seek medical attention if the patient experience chest pain, dizziness or shortness of breath or if you experience symptoms of overexertion.

WARNING:

Serious injury or death could result from operating the treadmill in the presence of explosive or flammable vapors and antiseptics.

Introduction

WARNING:

The potential for foot crush injury at frontal end of treadmill at lift mechanism (landing gear) when treadmill is descending.
Always Keep feet and hands away from this area.

Potential foot crush injury at rearward side rail, rear of side rail and rear roller exists when treadmill approaches full elevation.
Always Keep feet and hands away from this area.

Classification of Medical Device

This device is classified as follows, according to IEC 60601-1:

Medical Device Classification

Category	Classification
Type of protection against electrical shock	Class I motor operated physical medicine machine.
Degree of protection against electrical shocks	Type B external application applied part.
Degree of protection against harmful ingress or water	Ordinary equipment (enclosed equipment without protection against ingress of water).
Degree of safety of application in the presence of a flammable anesthetic mixture with air or with oxygen or with nitrous oxide	Equipment is not suitable for use in the presence of a flammable anesthetic mixture with air or with oxygen or with nitrous oxide.
Method(s) of sterilization or disinfection recommended by the manufacturer	Not applicable
Mode of operation	Continuous operation.

Regulatory and Safety Conformance

T2100-ST1 and T2100-ST2 treadmills meet FDA Class 1 motor operated physical medicine machines. They have been tested by Intertek Testing Services N.A Inc. and are listed by Engineering Testing Laboratories (ETL). However, the ultimate conformance to IEC 60601-1 is the responsibility of the system integrator when combined with other equipment. Additionally, all motorized equipment is potentially dangerous if used incorrectly. Before using the T2100-ST1 or T2100-ST2 treadmill, follow all precautions listed in this manual and read the entire Operator's Manual thoroughly. Use the T2100-ST1 and T2100-ST2 treadmills only as described.

NOTE:

This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when

Introduction

the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at owner's expense.

Table 1: Guidance and Manufacturer's Declaration - Emissions

The T2100-ST Series is intended for use in the electromagnetic environment specified below. The customer or user of the T2100-ST Series should ensure that it is used in such an environment.

Emissions Test	Compliance	Electromagnetic Environment - Guidance
RF Emissions CISPR 11	Group 1	The T2100-ST Series uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF Emissions CISPR 11	Class B	The T2100-ST Series is suitable for use in all establishments, including domestic, and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonics IEC 61000-3-2	Class A	
Flicker IEC 61000-3-3	Complies	

Introduction

Table 2: Guidance and Manufacturer's Declaration – Immunity All ME Equipment and ME Systems

The T2100-ST Series is intended for use in the electromagnetic environment specified below. The customer or user of the T2100-ST Series should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
ESD IEC 61000-4-2	±8kV Contact ±15kV Air	±8kV Contact ±15kV Air	Floors should be wood, concrete, or ceramic tile. If floors are synthetic, the r/h should be at least 30%
EFT IEC 61000-4-4	±2kV Mains ±1kV I/Os	±2kV Mains ±1kV I/Os	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1kV Differential ±2kV Common	±1kV Differential ±2kV Common	Mains power quality should be that of a typical commercial or hospital environment.
Voltage Dips/Dropout IEC 61000-4-11	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	>95% Dip for 0.5 Cycle 60% Dip for 5 Cycles 30% Dip for 25 Cycles >95% Dip for 5 Seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the T2100-ST Series requires continued operation during power mains interruptions, it is recommended that the T2100-ST Series be powered from an uninterruptible power supply or battery.
Power Frequency 50/60Hz Magnetic Field IEC 61000-4-8	30 A/m	30 A/m, 50/60Hz	Power frequency magnetic fields should be that of a typical commercial or hospital environment.

Introduction

Table 3: Guidance and Manufacturer's Declaration – Immunity ME Equipment and ME Systems that is NOT Life-supporting

The T2100-ST Series is intended for use in the electromagnetic environment specified below. The customer or user of the T2100-ST Series should ensure that it is used in such an environment.

Immunity Test	IEC 60601 Test Level	Compliance Level	Electromagnetic Environment - Guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	(V1)=3Vrms (3V, 6V at ISM + Amateur Frequencies)	Portable and mobile communications equipment should be separated from the T2100-ST Series by no less than the distances calculated/listed below: $D=(3.5/V1) (\text{Sqrt } P)$ 150kHz to 80MHz $D=(3.5/E1) (\text{Sqrt } P)$ 80 to 800 MHz $D=(7/E1) (\text{Sqrt } P)$ 800 MHz to 2.5 GHz where P is the max power in watts and D is the recommended separation distance in meters. Field strengths from fixed transmitters, as determined by an electromagnetic site survey, should be less than the compliance levels (V1 and E1). Interference may occur in the vicinity of equipment containing a transmitter.
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.7 GHz	(E1)=3V/m (3 V/m at 80 - 2,700MHz, AM Modulation 9-28V/m, 385 - 6,000MHz, FM or Digital Modulation)	

Table 4: Recommended Separation Distances between portable and mobile RF Communications equipment and the T2100-ST Series ME Equipment and ME Systems that is NOT Life-supporting

The T2100-ST Series is intended for use in the electromagnetic environment in which radiated disturbances are controlled. The customer or user of the T2100-ST Series can help prevent electromagnetic interference by maintaining a minimum distance between

Introduction

portable and mobile RF Communications Equipment and the T2100-ST Series as recommended below, according to the maximum output power of the communications equipment.

Max Output Power (Watts)	Separation (m) 150kHz to 80MHz $D=(3.5/V1)(\text{Sqrt } P)$	Separation (m) 80 to 800MHz $D=(3.5/E1)(\text{Sqrt } P)$	Separation (m) 800MHz to 2.5GHz $D=(7/E1)(\text{Sqrt } P)$
0.01	0.11667	0.11667	0.23333
0.1	0.36894	0.36894	0.73785
1	1.1667	1.1667	2.3333
10	3.6894	3.6894	7.3785
100	11.667	11.667	23.333

Responsibility of the Manufacturer

Full-Vision, Inc., is responsible for the effects of safety, reliability, and performance of the treadmill only if the following conditions are met:

- Assembly operations, extensions, readjustments, modifications, or repairs are carried out by GE Healthcare which is authorized for service and installation.
- The electrical installation of the relevant room complies with the requirements of the appropriate local, state, and other government regulations.
- The equipment is used in accordance with the instructions for use.

Introduction

Responsibility of the Customer

The customer is responsible for providing appropriate desks, chairs, electrical wall outlets, network connections, and analog phone lines, and for locating any of the system components described in this manual in compliance with all local, state, and national codes.

The customer is solely responsible for the training, instruction, supervision and safety of all users of the T2100-ST Series treadmill, and to use it as intended by the manufacturer. This device is intended to be used as a motion appliance to facilitate cardiac or VO₂ medical evaluation.


- Read the operator's manual before operating the T2100-ST Series Treadmill.
- Assist in off-loading the patient in the event of abnormal or unexpected operation of the treadmill.
- If the treadmill is not responding properly, stop the treadmill, assist in removing the patient off the running belt, unplug the treadmill power supply, and seek factory authorized repair before attempting to restart the treadmill.
- Never allow children or pets near the machine without qualified adult supervision.
- Note the location of stop and/or emergency stop controls and their operation before starting a test or workout.
- This device is not intended for use by persons with reduced physical, sensory or mental capabilities, or lack of experience and knowledge unless they have been given supervision or instruction concerning use of the appliance by a person responsible for their safety.
- Verify the Patient and Operator both know how to stop the machine in the event of malfunction or emergency.
- Patient should not wear loose fitting nylon material when exercising on this treadmill to avoid generating electro-static discharge.
- Never attempt to remove any article of clothing while the running belt is moving.
- All persons on and around the treadmill must wear enclosed, protective footwear. Shoelaces must be tight and not drape as to cause a trip or catch hazard.

Sandals, flip flops, slippers and the like are not considered enclosed, protective footwear.

- Walk in the center of the running belt. Contact with the side rail and the moving belt could cause injury.

Introduction

- The Patient must always wear the strap while operating the T2100-ST Series Treadmill.

- The Patient must always wear stop Tether lanyard wrist strap while operating the T2100-ST Series Treadmill.

- Place the treadmill on a hard, level and unobstructed surface. See Chapter 3, "Assembly and Setup".
- Check input power cord connection and location for hazardous pinch points before use.
- Check input communications cord connection (if equipped) for proper interface with all equipment.
- Keep all cords clear of patient to avoid trip hazards.
- Never attempt to remove the motor pan hood or do electrical repairs yourself. Repairs should only be done by a factory authorized repair provider.
- Always unplug the T2100-ST Series Treadmill when servicing, inspecting or cleaning the treadmill.
- Routinely inspect the treadmill for loose parts.
- Inspect handrails and ensure they will support the patient properly.

Introduction

- Always start the running belt at its slowest speed before starting the patient test.
- Do not step onto belt when it is moving.
- Always slow the running belt to its minimum speed before stopping.
- Keep hands, feet, and clothing away from any moving parts.
- Verify no one is near the elevation mechanism before operating. Never put any part of the body under any part of a running treadmill.
- Never drop or insert objects into any opening.
- Never drape garments, hook-up leads, or other equipment over the side rails or drop objects on the belt while the T2100-ST Series Treadmill is running.
- Do not allow moisture or oils to accumulate on equipment, creating a slip hazard.

Product and Package Information

This section describes the location of the labels used on your device and its packaging. It also describes the symbols used on the labels.

Symbols


The following symbols may appear on the device or its packaging. Familiarity with these symbols assists in the safe use and disposal of the equipment. For equipment symbols not shown, refer to the original equipment manufacturers (OEM) manuals.

Symbols are used to convey warnings, cautions, prohibitions, mandatory actions, or information. Any hazard symbols on your device or packaging with markings in color indicates there is certain danger and is a warning. Any hazard symbols on your device or packaging that is in black and white indicates a potential hazard and is a caution.

Symbols

Symbol	Description
 ABC123	Catalog or Orderable Part Number Indicates the manufacturer's catalog or part number.
 ABC123	Serial Number Indicates the manufacturer's serial number.

Introduction

Symbol	Description
MANUFACTURED BY:  YYYY-MM-DD	Manufacturer Name, Address and Manufacturing (Year-Month-Date) Indicates the name and address for the manufacturer of this device.
	CAUTION: CONSULT ACCOMPANYING DOCUMENTS - There may be specific warnings or precautions associated with the device that are not otherwise found on the label. Consult the accompanying documentation for more information about safely using this device.
	CAUTION: ELECTRIC SHOCK - Indicates the presence of hazardous energy circuits or electric shock hazards. To reduce the risk of electric shock hazards, do not open this enclosure. Refer servicing to qualified personnel.
	Reading of the Owner's Manual is mandatory.
	Upper Temperature Limit Indicates the maximum temperature for transportation and handling of this package.
	Temperature Limits Indicates the upper and lower temperature limitations for the transportation and handling of this package.
	Humidity Limits Indicates the upper and lower Non-Condensing Humidity limitations for the transportation and handling of this package.

Introduction

Symbol	Description
	European Union Disposal Requirements This equipment complies with the EU WEEE marking requirement for proper disposal of electrical and electronic waste in accordance with the European Directive. This directive calls for separation and recovery or reuse of used electrical or electronic equipment upon end-of-life EEE disposal. The T2100 must not be disposed of as unsorted municipal waste. Electrical or electronic components must be collected separately and disposed of in accordance with your local requirements and sources. The EEE program minimizes any potential effects on the environment and user health by eliminating the potential presence of hazardous substances in the waste stream. Customers should contact their local authorities or T2100 Distributor for guidance in complying with the directive.
	Keep Dry Indicates that you need to keep the container away from rain and other sources of moisture.
	CE Mark Indicates the device or product conforms with applicable EU (European Union) directives.
	Eurasian Conformity mark Conformity to applicable technical regulations of Customs Union.



Introduction

Symbol	Description
	Electrical Testing Laboratories Indicates the device or product has been tested by an accredited third-party testing laboratory and meets applicable safety standards for sale and distribution within North America.
	PCT (GOST-R) Mark Indicates the device or product conforms with applicable Russian Ghost standard technical and safety standards.
	Protective earth (ground).
	Alternating current.
	Device is suitable for the external application of the type "B" applied parts.
	Unique Device Identification is a unique marking of the medical device
	Follow operating instructions
	Medical device
	Authorized Representative in European Union (Regulatory affairs only).
	Authorized Representative In Switzerland
	United Kingdom Responsible Person
	Federal law restricts this device to sale by or on the order of a physician.
	Importer: Indicates the entity importing the medical device into the locale.

Label Locations


This section identifies the labels and their locations on the product and packaging.

Introduction



DETAIL 1

Introduction



Detail 2




Refer to the previous illustrations for the locations of the labels identified in the following table. For detailed descriptions of the symbols that appear on the labels, refer to "Symbols" on page 22.

Introduction


Item	Label	Location	Description
1		TMX428 on center handrail TMX428CP center console	Identifies Reading of the Owner's Manual
2		Stop on RH or LH end cap on handrail	Identifies E-Stop
3		Front of the device on the hood	Identifies the Caution Electrical shock hazard
4		On the motor	Identifies the Caution Electrical shock hazard
5		On the Power Cord	Identifies DC Hi-Pot Caution
6		Front of the device on the hood	Identifies Master Power switch.
7		Front of the device	Identifies Power Off (disconnection form mains) and Power On (connection to the mains)
8		Front of the device	Identifies the location for the RS232 & USB connection point.
9	Decal T2100-STXN	On side channel Right & Left position above rear foot.	Identifies GE Model T2100-STX
10		Front of device	Identifies the product model.
11		Front of device	Identifies Listing Standards
12		Front of device	Contains the European Union disposal requirements

Introduction

13	 001 00860176000637 111 190101 211 FVDC-XXXX	Front of the device on motor pan & on Shipping label.	Identifies Unique Device Identifier (01) GTIN Number (11) Manufacture Date (21) Serial Number
14	 T2100-STX 2097357-XXX GDC-XXXX 4TH EDITION PRODUCTION 2097357-XXX REV XX FULL VISION 2097357-XXX	On shipping container.	Identifies the following information for shipping: Customer Model Number Revision Serial Number Storage conditions Regulatory Compliance Country of origin EC Representative information
14A	 INFORMAÇÕES DO IMPORTADOR E DISTRIBUIDOR INFORMAÇÕES DO FABRICANTE INFORMAÇÕES DO PRODUTO Segurança Computadorio	On shipping container T2100-ST2 Only.	Identifies Brazil shipping registration information
15		Front of device	Identifies Medical Device
16		Front of device	Identifies FSE Barcode for service information requested by GE Healthcare
17		Front of device	Identifies 4TH Edition Production and Revision level
18		Front of device	Identifies Japan registration information

19		Front of device & On shipping container T2100-ST2 Only.	Identifies Indonesia registration information
20		Front of device	Identifies Malaysia registration information
21		Front of device	Identifies Brazil registration information
22		Front of device	Identifies China registration information
23		Front of device	Identifies Authorised REP information
24	 	On shipping container T2100-ST2 Only.	Identifies Indonesia & India registration information

Product Label

Item	Name	Description
1	Manufacturer	Full Vision Inc.

Introduction

Item	Name	Description
2	Private Label Manufacturer	GE Medical Systems Information Technologies, Inc
3	Model Number	Identifies model of treadmill
4	Manufacturer Model Number	Identifies manufacturing model of treadmill
5	Serial Number	Manufacturer assigned serial number
6	REF	GE Medical Systems reference part number
7	Manufacturer Part Number	Manufacturing part number
8	Manufacturer Date	Manufacturing date code
9	Voltage	Specifies operating voltage of treadmill
10	Hertz	Specifies the electrical hertz of treadmill
11	Amps	Specifies amperage of treadmill
12	Type B Equipment	Device is suitable for the external application of type "B" applied parts

Service Information

This section provides information pertaining to the maintenance and servicing of the system. Familiarize yourself with this information before requesting service from GE Healthcare or its authorized representatives.

Service Requirements

Failure on the part of the responsible individual, hospital, or institution using this equipment to implement a satisfactory maintenance schedule may cause undue equipment failure and possible safety hazards.

Regular maintenance, irrespective of usage, is essential to ensure that the components of this system are always functional when required.

Warranty Information

This device is considered GE Healthcare-supplied hardware. Only authorized GE Healthcare service personnel should service the device. Any unauthorized attempt to repair equipment under warranty voids that warranty. It is the user's responsibility to report the need for service to GE Healthcare or to one of their authorized agents.

Additional Assistance

GE Healthcare maintains a trained staff of application and technical experts to answer questions and respond to issues and problems that may arise during the installation, maintenance, and use of this system.

Introduction

Contact your local GE Healthcare representative to request additional assistance.

Reporting of serious incidents

Any serious incident related to the use of this GE Healthcare device, should be reported to both, the manufacturer and the health authority/competent authority where the device is installed.

- To report to manufacturer:

Email: tmsservice@full-vision.com

- Please provide the following information:
 - The catalogue number or the model designation of the device as stated on its identification plate affixed on the device
 - The SystemID/serial number/lot number of the device
 - Date of incident
 - Description of incident, including any patient or user impact/injury
 - Your contact information (facility, address, contact name, title, and telephone number)

Manual Information

This section provides information for the correct use of this manual.

Keep this manual with the equipment at all times and periodically review it. You should request training assistance from GE Healthcare, if needed.

Manual Purpose

This manual provides information necessary for the configuration and safe operation of this equipment in accordance with its function and intended use. It is not intended as a replacement for, but a supplement to, thorough product training. Keep it with the equipment at all times. Additional manuals may be ordered by contacting GE Healthcare.

Refer to the service manual for technical information related to the maintenance and repair of the equipment.

Document Conventions

This document uses the following conventions.

Introduction

Typographical Conventions

The following table identifies the typographical conventions used in both this document and GE Healthcare Diagnostic Cardiology product documents.

Convention	Description
Bold Text	Indicates keys on the keyboard, text to enter, or hardware items such as buttons or switches on the equipment
<i>Italicized Bold Text</i>	Indicates software terms that identify menu items, buttons, or options in various windows.
KEY1+KEY2	Indicates a keyboard operation. A plus (+) sign between the names of two keys indicates that while holding the first key, you should press and release the second key. For example, Press CTRL+ESC means to press and hold the CTRL key and then press and release the ESC key.
<space>	Indicates that you must press the spacebar. When instructions are given for typing a precise text string with one or more spaces, the point where you must press the spacebar is indicated as: <space>. This ensures that the correct number of spaces are inserted in the correct positions within the literal text string. The purpose of the <> brackets is to distinguish the command from the literal text within the string.
Enter	Indicates that you must press the Enter or Return key on the keyboard. Do not type Enter .
>	The greater than symbol, or right-angle bracket, is a concise method to indicate a sequence of menu selections. For example, the statement "From the main menu, select System > Setup > Options to open the Option Activation window" replaces the following: 1. From the main menu, select System to open the System menu. 2. From the System menu, select Setup to open the Setup menu. 3. From the Setup menu, select Options to open the Option Activation window.

Illustrations

All illustrations in the document are provided as examples only.

Notes

Notes provide tips or additional information that, while useful, are not essential to the correct operation of the tools. They are called out from the body text through a flag word and indentation, as follows:

NOTE:

The tip or additional information appears indented below the **NOTE** flag word.

Introduction

Related Documents

The following documents are referenced in this manual and provide additional information that may be helpful in the installation, configuration, maintenance, and use of this product.

Part Number	Title
2097937-001	T2100-ST1 /T2100-ST2 Operator's Manual, English
2097937-021	T2100-ST1/T2100-ST2 Operator Manual, S. Chinese
2097937-051	T2100-ST1/T2100-ST2 Operator Manual, Danish
2097937-061	T2100-ST1/T2100-ST2 Operator Manual, Dutch
2097937-091	T2100-ST1/T2100-ST2 Operator Manual, French
2097937-131	T2100-ST1/T2100-ST2 Operator Manual, Indonesian
2097937-151	T2100-ST1/T2100-ST2 Operator Manual, Japanese
2097937-171	T2100-ST1/T2100-ST2 Operator Manual, Korean
2097937-211	T2100-ST1/T2100-ST2 Operator Manual, Portuguese Brazil
2097937-271	T2100-ST1/T2100-ST2 Operator Manual, Spanish
2097937-011	T2100-ST1/T2100-ST2 Operator Manual, Bulgarian
2097937-031	T2100-ST1/T2100-ST2 Operator Manual, Croatian
2097937-041	T2100-ST1/T2100-ST2 Operator Manual, Czech
2097937-071	T2100-ST1/T2100-ST2 Operator Manual, Estonian
2097937-081	T2100-ST1/T2100-ST2 Operator Manual, Finnish
2097937-101	T2100-ST1/T2100-ST2 Operator Manual, German
2097937-111	T2100-ST1/T2100-ST2 Operator Manual, Greek
2097937-121	T2100-ST1/T2100-ST2 Operator Manual, Hungarian
2097937-141	T2100-ST1/T2100-ST2 Operator Manual, Italian
2097937-161	T2100-ST1/T2100-ST2 Operator Manual, Norwegian
2097937-181	T2100-ST1/T2100-ST2 Operator Manual, Polish
2097937-201	T2100-ST1/T2100-ST2 Operator Manual, Portuguese EU
2097937-191	T2100-ST1/T2100-ST2 Operator Manual, Romanian
2097937-221	T2100-ST1/T2100-ST2 Operator Manual, Russian
2097937-231	T2100-ST1/T2100-ST2 Operator Manual, Serbian
2097937-241	T2100-ST1/T2100-ST2 Operator Manual, Slovak
2097937-261	T2100-ST1/T2100-ST2 Operator Manual, Swedish
2097937-251	T2100-ST1/T2100-ST2 Operator Manual, Turkish
2097937-281	T2100-ST1/T2100-ST2 Operator Manual, Vietnamese
2097937-311	T2100-ST1/T2100-ST2 Operator Manual, Lithuanian
2097937-291	T2100-ST1/T2100-ST2 Operator Manual, Kazakh
2097937-301	T2100-ST1/T2100-ST2 Operator Manual, Latvian

Introduction

Part Number	Title
2097937-321	T2100-ST1/T2100-ST2 Operator Manual, Slovene
2097937-331	T2100-ST1/T2100-ST2 Operator Manual, Ukrainian

Introduction

Training

This manual is intended as a supplement to, not a substitute for, thorough product training. If you have not received training on the use of the system, you should request training assistance from GE Healthcare.

To see available training, go to the GE Healthcare training Web site (<http://www.gehealthcare.com/user/education/index.html>). Under the Technical Service Education section, select Diagnostic Cardiology.


For more self-paced course offerings, tools, and reference guides you may find useful, please visit the GE Healthcare Education Store at www.gehealthcare.com/educationstore.

2

Product Overview

The T2100-ST1 and T2100-ST2 treadmills are designed and built to withstand the extraordinary demands of medical devices and are compatible with CASE, CardioSoft/CS, MAC 5500, and MAC 2000.

References to left, right, front, and rear assume that you are standing on the treadmill, facing the handrails. All parts listed below are considered Patient Applied Parts except where noted.

Item	Description
1	Emergency Stop Button
2	Patient Handrails
3	Pull Tether Switch
4	Side Rail
5	Running Belt
6	Elevation Landing Gear (Non-Applied Part)
7	Rear Foot (Non-Applied Part)

Product Overview

Safety Systems

- Dual comparative speed sensors
- Auto runaway shutdown
- Auto communication loss shutdown
- Manual twist lock Emergency Stop button
- Manual Stop Tether
- Braking system for safe patient off-loading
- Fire-rated motor pan hood enclosure

Treadmill

- Patient weight capacity 500 lb., 227 kg
- All steel construction with powder-coat finish
- Treadmill net weight: 425 lb., 193 kg

Drive System

- Heavy-duty 4-peak hp. brushless, DC servo motor
- T2100-ST1 110-120VAC, single-phase, 60 Hz, 20-amp power supply
- T2100-ST2 200-240VAC, single-phase or split phase, 50-60 Hz, 15-amp power supply

Speed Range

- 0.1 to 15.0 mph, 0.2 to 24.0 km/h, self-calibrating and adjustable in 0.1 mph 0.1 km/h increments.

NOTE:

The T2100-ST2 maximum speed (15.0 mph/24.0 km/h) will deteriorate at lower voltages (210VAC or below).

Incline Range

- 0 to 25%, 0.5% incremental movements, self-calibrating.

Running Surface

- 22in. x 63in. / 56cm x 160cm
- MasterTrack® running belt tracking system
- Cushioned running deck absorbs shock of foot falls
- Self-lubricated and reversible running deck
- Step-up height (7 in. / 18cm from floor)

Product Overview

Communication Ports

- RS232 Female Serial port
- USB 1.0 "B" port

Floor Surface Footprint

- 33.0 in x 78.5 in, 84cm x 200cm level surface. (See "Location" on page 48.)

Operating and Storage Condition Recommendations

- Operating Temperature Range: 4.5° to +38° C (+40° to +85° F)
- Storage Temperature Range: -40° to +70° C (-40° to +158° F)
- Operating and Storage Relative Humidity Range: 10% - 90%, non-condensing
- Altitude: -50 to 9,842 feet (-15m to 3000m); derate performance by 5% per each additional 500 feet (152m) above 5,280 feet (152m)

Power Requirements

The T2100-ST1 is designed to operate on dedicated 100-130 VAC 20-amp power, and the T2100-ST2 is designed to operate on dedicated 200 to 264 VAC 13-amp power. Make sure that the treadmill is connected to an outlet that looks like one of the following illustrations.

- T2100-ST1 Power Consumption: 1934 watts (6600 BTU), 20 amps
- T2100-ST2 Power Consumption: 2580 watts (8804 BTU), 13 amps

This product is equipped with a three-wire grounding-type plug. The plug will only fit into a grounding-type outlet. This safety feature must not be disabled. Contact a qualified electrician if you are unable to insert the plug into your outlet, or uncertain if the outlet meets local electrical codes. Polarized outlets such as NEMA 5-20, CEE7/7, and BS1363 must be verified for proper polarity configuration before plugging in the T2100-ST1 or T2100-ST2. Incorrect polarization of the outlet could cause failure of onboard electrical components or cause electrical shock. Proper grounding is necessary for the equipment to meet acceptable current leakage standards consistent with the standards to which it was certified.

Product Overview

110-120 VAC
NEMA 5-20R
Single Phase


220/240 VAC
NEMA 6-15R
Split Phase


230 VAC
AS/NZS 3112


220/240 VAC
CEE 7/7 EURO


250 VAC
Type N BRAZIL


230 VAC
BS1363


240 VAC
BS546 3 PIN

250 VAC
Type K DANISH

240 VAC
GB 1002


250 VAC
Type H ISREAL

250 VAC
Type L ITALY

250 VAC
Type J SWISS

250VAC
BS546 INDIAN

250 VAC
NEMA 6-20P
Single Phase

WARNING:

The T2100-ST1 and T2100-ST2 treadmills must be grounded to reduce the risk of electrical shock. If a malfunction occurs, earth grounding provides a path of least resistance for an electric current. Ungrounded connections must not be used.

No other equipment may be used on the electrical circuit with the treadmills. Do not use extension cords. Using a shared or unreliable circuit can also cause the treadmills to

Product Overview

unexpectedly shut off, potentially resulting in injury to the patient.

Ensure the power switch is in the off position before plugging in the T2100-ST1 or T2100-ST2. A power surge could damage the sophisticated electronic system of the treadmills.

NOTE:

The T2100-ST1 and T2100-ST2 Treadmills must have their own dedicated power outlet.

Maximum length of power cord should not exceed 10 Feet or 3.05 Meters. Replacement Power Cords need to be ordered through the manufacturer of the operating equipment.

WARNING:

Use of this equipment adjacent to or stacked with other equipment should be avoided because it could result in improper operation. If such use is necessary, this equipment and the other equipment should be observed to verify that they are operating normally.

WARNING:

Only use accessories, transducers, and cables provided by the manufacturer of the equipment. Using alternative non supplies accessories of the components listed above could result in increased electromagnetic emissions or decreased electromagnetic immunity of the equipment and could result in improper operation.

3

Assembly and Setup

The T2100-ST1 and T2100-ST2 treadmills are shipped fully assembled and packaged in a folded condition. They are designed to pass through a standard 36" door opening measuring at least 35½". It will be necessary to remove the door from the jam in most cases if the door is not capable of opening fully parallel to door opening. After you have unpacked the treadmill and secured the handrail assembly to the frame, move the treadmill to the area by rolling it on its front wheels. If your treadmill must pass through a door opening less than 36" wide, additional disassembly by removal of the handrails is required. Refer to "Removing and Reinstalling the Handrails for Moving" on page 190. This task should be performed by an authorized service provider to ensure that the T2100-ST1 and T2100-ST2 are properly reassembled and functioning correctly.

NOTE

The treadmill arrives already calibrated for speed and elevation in US Voltage configuration. The speed and elevation may need to be calibrated to accommodate local installation voltage.

Safe Handling Guidelines

WARNING

The T2100-ST1 and T2100-ST2 treadmills each weigh 425 lbs. When lifting the rear of the treadmill and rolling it on its front

Assembly and Setup

wheels, you are moving 132 lbs. This requires 2 people to safely remove it from the pallet and locate it in the facility.

If you are moving the treadmill over rough surface, such as pavement, use a dolly under front of the treadmill to prevent damage to the wheels and lift mechanism.


- Do not attempt to move the treadmill with the handrails in the shipping position due to the possibility of cutting the internal wiring. You must either fully secure the handrails in their proper position or secure handrails with 3/8-16 bolts and 3/8 lock washers in the folded position.
- Lift the end of the bed assembly to a comfortable height, keeping knees bent and back straight as you lift.
- Rotate the treadmill in the direction you want to go (the treadmill will pivot on its wheels) and push forward.
- When you have maneuvered the treadmill into its location, gently lower the end of the bed assembly to the floor.

Uncrating Instructions

Tools required for uncrating:


- Safety Knife

1. The treadmill is shipped in a cardboard box with a wood pallet. Cut the 3 nylon straps holding the cardboard cover to the wood pallet.

2. Remove the cardboard box by lifting straight up from one end to clear the treadmill.


Assembly and Setup

3. Carefully cut the zip tie to open the protective plastic bag.
4. Pull the plastic bag down and fold over all 4 corners of the pallet.

Note: This will prevent the plastic bag being rolled up into the elevation wheels when removing from wood pallet.

5. Carefully cut the nylon strap securing the treadmill to the pallet.

6. DO NOT REMOVE nylon strap for the handrails.

CAUTION

DO NOT REMOVE treadmill from wood pallet with a forklift. The treadmill may slide off fork causing damage to handrails & side channels.


WARNING

The device weight is 425lbs. This requires 2 people to safely remove it from the pallet and locate it in the facility.

7. Remove from the wood pallet by lifting from the rear of the treadmill and slide off one corner at a time.

Assembly and Setup

Step 3


NOTICE


- When moving leave handrails in folded position secured with nylon strap.
 - **NOTE** ▪ Required minimum clearance of 35 1/2" (902mm) to prevent handrail damage.
- At final installation set-up handrails and remove packaging materials.
 - **NOTE** ▪ Removing packaging materials with sharp objects may cause cosmetic damage to treadmill.

8. When moving or relocating the treadmill use a furniture dolly or floor dolly.
9. Position the dolly at the elevation pivot point to provide treadmill balance.

Furniture Dolly or Floor Dolly


Note: When moving to final destination you may experience a narrow passage or doorway. The treadmill can be tipped on the side with handrails in the folded position. This will allow clearance to a minimum of 28" (711mm) to safely pass thru opening.

10. Remove nylon strap allowing the handrail to be installed in upright position

Assembly and Setup

(See Initial Setup Instructions below).

11. Carefully remove the stretch wrap and protective packaging material starting at the top of the handrail.
12. Dispose of packaging material by approved local requirements.


Initial Setup

Tools required for assembly:

- 5/16 Allen wrench (supplied)

The treadmill is shipped with the handrails loose, straddling the treadmill frame. It is advised that you secure the handrails in their proper location before removing the treadmill from the base of the crate. This prevents the internal wires running down the handrail mount to the motor pan from being cut.

1. Swing the handrail assembly into the operating position and insert (2) 3/8-16 bolts and 3/8 lock washers on each side and tighten securely.
2. Install (2) plastic caps each side for a finished look.


Assembly and Setup

Item	Description
1	Operating Position
2	Shipping Position
3	Insert (2) washer and bolts each side
4	Insert (2) caps each side
5	Pivot Point
6	When folding handrails apply cardboard between frame and handrail to prevent handrail damage.

Location

Place the T2100-ST1 or T2100-ST2 on a firm and level hard surface that is free of tile grout lines. The illustration below shows the minimum recommended clearances from the treadmill edges to any obstruction for dismount and safety purposes. Observe that the operator should be stationed by the Emergency Stop Button (ESB).

WARNING:

The T2100-ST1 and T2100-ST2 treadmills conform to FCC class B rating for electromagnetic emissions. It is recommended not to place the treadmill closer than 5ft. (1.5m) from sensitive electronic devices within the room or in an adjacent room. If an interference problem occurs, move the treadmill farther away from the sensitive device or relocate either device to another area, or consult with an EMI specialist for ways to shield the room from electromagnetic radiation.

Do not place it on thick or long-pile carpeting. Such carpeting could cause instability or static build-up, and carpet fibers could get caught in the belt and damage the unit.

Assembly and Setup

Ensure that power cords do not cross traffic areas. Exposed power cords can cause a fall, resulting in injury.

Keep it away from sources of moisture, such as spas or fountains. Moisture can cause the electronic circuitry to malfunction.

Electrical Safety Tests

Except where noted, the following electrical safety tests are optional at the time of install.

You must conduct the following electrical safety tests after any service is performed or if the hood was removed.

AC Line Voltage Test

This test verifies that the wall outlet supplying power to the equipment is properly wired. For international wiring test, refer to the internal standards agencies of that country.


NOTE

The AC Line Voltage test is required at the time of installation.

100 to 130 VAC, 60 Hz

Using a digital voltmeter set to measure at least 200 VAC to check the voltage of the NEMA 5-20R AC wall outlet (U.S. only or applicable international connection; dedicated 20-amp service). If the measurements are out of range, have a qualified electrician repair the outlet. The voltage measurements should be:

100 to 130 nominal VAC between the "neutral" and "hot" contacts.


Assembly and Setup

200 to 264 VAC, 50/60 Hz

Using a digital voltmeter set to measure at least 300 VAC to check the voltage of the NEMA 6-15R AC wall outlet (U.S. only or applicable international connection; dedicated 13-amp service). If the measurements are out of range, have a qualified electrician repair the outlet. The voltage measurements should be:

200 to 264 nominal VAC between the two "hot" contacts.

Power Cord and Plug

Verify the power cord being used with the treadmill is good:


- Verify that the line, neutral (if applicable), and ground conductors are properly connected to the power cord plug and are not short-circuited. Replace the power cord if necessary.
- Failure of the power cord strain relief is very common. Often users of the equipment pull on the power cord itself, rather than the power cord plug, to unplug the unit from a wall receptacle. If in doubt, test for continuity through each.
- Perform the Ground Continuity Test, or the test method that is required by your Country/Local governing safety organization. For international power cords, refer to the internal standards agencies of that particular country.

Assembly and Setup

Ground Continuity Test

This test verifies that there is continuity (less than 100 mΩ resistance) between exposed metal surfaces, which have the potential to become energized, and the ground prong on the mains AC power cord. Look for an exposed metal screw, or, if the metal surfaces are anodized or painted, scrape off a small area in an inconspicuous area on the aluminum casting, for the probe to make direct contact with the metal.

1. Disconnect the power cord from wall outlet with the T2100-ST Series treadmill power switch in the "ON" position.
2. Use a digital multimeter to check metal surfaces of the equipment as illustrated below. Make adjustments for any resistance in the test leads.

If the measurements are significantly out of range, check for breaks in the power cord or in the internal connections within the units.

Conducting Leakage Tests

The leakage tests are safety tests to ensure that the equipment poses no electrical health hazards. Use the table below to determine which tests apply to the unit under test and the maximum allowable leakage currents. For international leakage limits, refer to the internal standards agencies of that country.

If the unit under test fails the leakage tests, do not allow the customer to use the equipment. Call Tech Support for assistance.

GE Healthcare recommends that you perform these tests:

- Before applying power for the first time
- Whenever internal assemblies are serviced

Assembly and Setup

NOTE:

The accuracy of the leakage tests depends on a properly wired wall outlet. Do not proceed until you verify the integrity of the power source.

Earth/Chassis Leakage Current to Ground					limits
1.	Forward Polarity (L1, L2)	NC	_____ μA	Pass/Fail	500 μA
2.	Neutral (L2) open, Forward Neutral Polarity	SFC	_____ μA	Pass/Fail	500 μA
3.	Ground open, Forward Polarity	SFC	_____ μA	Pass/Fail	500 μA
4.	Ground open, Reverse Polarity	SFC	_____ μA	Pass/Fail	500 μA
5.	Neutral (L2) open, Reverse Polarity	SFC	_____ μA	Pass/Fail	500 μA
6.	Reverse Polarity (L2, L1)	NC	_____ μA	Pass/Fail	500 μA
Ground Continuity					Resistance
1.	AC mains power cord ground prong to exposed metal surface (ground lug)	N/A	_____ $\text{m}\Omega$	Pass/Fail	100 $\text{m}\Omega$

WARNING:

Total system leakage current must not exceed 500 microamperes.


Leakage Test Diagrams

These diagrams show only a representation of how a typical leakage current tester function. Follow the instructions provided with the leakage current tester that you use.

Test #1

Ground-Wire-Leakage-to Ground

Make sure the UUT power switch is in the ON state without running belt running.


Assembly and Setup

Test #2

Chassis-Leakage-to-Ground (Exposed Chassis)

Make sure the UUT is in the ON state with running belt running at 2.0mph.

Secure Cables

Tie down cables to ensure they do not get caught in the wheels or the elevation racks.

Running Belt Tracking Adjustment

NOTE:

This adjustment is not covered under your warranty after installation. It is important that you review these instructions thoroughly before proceeding with belt tracking adjustment. Uneven floors accelerate belt misalignment. This situation may require more frequent adjustment to prevent belt damage.

The patented MasterTrack® Belt Tracking System significantly reduces the need to adjust the belt on your T2100-ST1 or T2100-ST2 treadmill. However, when you operate your treadmill for the first time, you may need to adjust the tracking of the belt to conform to your floor. You may also need to adjust the tracking if you move the machine to another location. (See "Running Belt Tracking Adjustment" on page 71 for details.)

Running Belt Tension Adjustment

Your running belt has been pre-tensioned at the factory and run for 16 hours prior to shipment. It may, however, be necessary to adjust the belt tension when the treadmill is run in its final location. A loose belt tends to hesitate or stick with a heavy foot plant. If your belt needs tensioning, the adjustment procedure can be found in "Running Belt Tension Adjustment" on page 72.

NOTE:

Improper adjustment could cause the treadmill to hesitate and cause a trip and fall hazard. This adjustment is not

Assembly and Setup


covered under your warranty after installation. It is important that you review these instructions thoroughly before proceeding.

Drive Belt Tension Adjustment

The drive belt tension has been pre-set at the factory to minimize maintenance. If there are indications that the drive belt has stretched and become loose, refer to Drive Belt Tension Adjustment " on page 73. Symptoms of a stretched drive belt could include increased noise.

Test Plug Procedure

Each T2100-ST1 and T2100-ST2 treadmill includes an RS-232 test plug that enables you to test the operation of the treadmill without the host controlling device attached. The plug is located on the left side of the treadmill secured to the frame by Velcro®. The plug is to be used only for testing the treadmill. Do not stand on or use the treadmill while testing.


Item	Description
1	Test Plug Location on Left
2	Circuit Breaker Array
3	Incoming Power Cord
4	Main Power Switch

Assembly and Setup

Communication Ports

The communication ports are located at the very front of the treadmill near the center of the unit.

Two ports are offered with equal communication capability. The standard female RS232 port and a USB Type B port offer connectivity diversification.

Item	Description
1	RS232 Port "Female"
2	USB "B" Port

Using the test plug

1. Turn the power "OFF" at the treadmill.
2. If connected, remove the RS232 or USB interface cable from the treadmill and plug in the test connector into the RS232 port.
3. Press and hold the button on the test connector and turn the treadmill power "ON". Continue holding until the treadmill begins to elevate.
4. When the treadmill begins to raise, repeatedly press the button to cycle through the different phases of the test.

Assembly and Setup

The following table shows the sequence of self-test stages:

Stage	Incline	Speed (mph)
1	5%	0.0
2	10%	0.0
3	15%	0.0
4	20%	0.0
5	20%	1.0
6	20%	2.5
7	20%	5.0
8	20%	7.5
9	20%	10.0
10	15%	7.5
11	10%	5.0
12	5%	2.5
13	0%	1.0
14	0%	0.0

NOTE:

Successful completion of the preceding testing procedure ensures that the treadmill is fully functional and responsive to command signals.

NOTE:

Unsuccessful completion of the preceding testing procedure indicates a problem with the setup. Call GE Medical Systems Information Technologies to trouble shoot failure of test plug procedure.

5. Turn the power "OFF" at the treadmill.
6. Remove the test connector.
7. Reconnect the RS232 or USB interface cable from the host computer.

You are ready to begin the set-up procedure prescribed by your medical test equipment supplier.

Assembly and Setup

Setup and Connection to Host

Currently, four GE host stress systems can connect to the T2100-ST series treadmill and are used as controlling devices: CASE, CardioSoft/CS,:

Note: The treadmill with serial number GEDC-6608 onwards (firmware version 1.19.109) are not compatible with MAC 5500ST and MAC 2000ST

CASE connection and setup:

1. Connect one end of the RS232 serial cable (GE part number 2097829-109) to the RS232 COM port located on the lower front panel of the T2100-ST series treadmill and run it through the communication cable clamp, as shown in the following photograph.

2. Connect the other end of the RS232 serial cable (GE part number 2097829-109) to a COM port on the CASE unit (it is suggested to use COM1 or COM2).
Referencing the CASE rear I/O connector panel label, take note of the COM port number plugged into.
3. Make sure the treadmill is connected to AC power and turned on.
4. If necessary, turn on CASE with the power button on the side panel.
The program will start up automatically and the Windows screen opens.
5. Log in when prompted.
Refer to the CASE Operator's Manual for logon instructions.
6. If not already in the software application, double-click the CASE icon.
The CASE initial screen opens.

Assembly and Setup


CASE Initial Screen

NOTE:

Depending on the version of software and user preferences, the initial screen will vary.

7. On the upper right side of the screen, click the **System Configuration** button.

Depending on the version of software the button may look like one of the following:

The **System Configuration** window opens.

8. Click the **Devices** tab.

The COM port number is displayed in the **Port** field, as shown in the following screen.

Assembly and Setup

System Configuration

General | Devices | **Modem** | MUSE | Option Code | Country Settings

Treadmill Model	T2100	Ergometer Model	Internal
Port	COM 1	Port	COM A
BP Monitor	no	SpO2 Monitor	no
Port	COM 1	Port	COM 1
☐ Metabolic Cart		Ambulatory Blood Pressure	no
	COM 2	Ambulatory BP Device Port	COM 2
		Spirometry Device Type	no
		Spirometry Device Port	COM 2

Path to Hoho ECG Program

Path to External Program

☐ Use External Program for Spirometry

Procedure Name for External Program

Analog 1	none	TTL output	
Analog 2	none	Pulse Width (ms)	50
Analog 3	none	Delay (ms)	0
Analog 4	none	Polarity	positive

Print Help OK Cancel

9. Ensure the **Treadmill Model** is set to **T2100**.
10. Ensure the COM port number is set to the one physically connected to in step 2, above, and, if necessary, change to match.
11. Click **OK** to close the **System Configuration** window and return to the CASE initial screen.

If needed, the CASE device operating instructions should be referenced for controlling the treadmill with the CASE application.

CardioSoft/CS connection and setup:

To communicate to the T2100-ST1 or T2100-ST2 treadmill using the USB port, you will need to install the appropriate USB driver software on your host computer. The USB driver is supplied on the provided flash drive. When connecting to the USB port, ensure port configuration is congruent with your software port identification.

First Time Installation of the Serial Driver Software on the CardioSoft/CS PC

This section installs the USB Serial Driver on the PC which has the CardioSoft or CS application installed.

NOTE:

Do not plug any hardware, including the treadmill USB cable, into the CardioSoft/CS PC until directed to do so.

Assembly and Setup

1. If necessary, turn on the CardioSoft/CS PC and log on when prompted.
Refer to the *CardioSoft/CS Operator Manual* for logon instructions.
2. Locate the flash drive (PN 317-160-313) shipped with the T2100-ST series treadmill and insert into a USB port of the CardioSoft/CS PC.
3. Browse the flash drive and locate the file named **CDM21218_Setup.exe** and double-click it to start the installation.

NOTE:


Do not download the driver from the FTDI website.

NOTE:

The following instructions are for a PC running Windows 7. The basic instruction flow should be the same for Windows 8 and 10, although the screens may vary. If necessary, refer to Appendix 0, "Reference Documents", for a list of reference documents that contain helpful OS-specific guidance.

4. If you receive the prompt "**Do you want to allow the following program to make changes to this computer?**", verify the program name, **CDM21218_Setup.exe**, and publisher, **Future Technology Devices International Ltd**, and click **Yes**.

The **Driver Extract** window opens.

5. Click **Extract**.

Assembly and Setup

The Device Driver Installation Wizard opens.

6. Click **Next**.

The **License Agreement** window opens.

7. Select **I accept this agreement** and click **Next**.

An installation progress window opens.

Assembly and Setup

The installation process should take less than a minute. When the installation is complete, the **Device Driver Installation Wizard** window opens confirming a successful installation.

8. To close the window, click **Finish**.

Assembly and Setup

Communication Port Identification and CardioSoft/CS setup

1. Connect a USB cable (GE part number 2044095-001) between the CardioSoft/CS PC and the T2100-ST series treadmill.
2. Connect AC power to the treadmill.
3. Power on the treadmill.
4. If necessary, turn on the CardioSoft/CS PC and the monitor.
5. On the CardioSoft/CS PC, determine the COM port that the USB Serial Port is using by doing the following:

- a. On the CardioSoft/CS PC click **Start** and in the **Search programs and files** field, type: **Device Manager** and press **Enter**.

The **Device Manager** window opens.

- b. Scroll to **Ports (COM & LPT)**.
- c. Find the **USB Serial Port** in the list

The port number is given in the parenthesis.

NOTE:

This may take a minute or so to display.

- d. Make note of this port and port number, you will use it later when selecting the COM port in the **System Configuration** window of the CardioSoft/CS application software.
- e. Close **Device Manager** window.
6. If not already in the software application, double-click the CardioSoft or CS icon.

The CardioSoft or CS initial screen opens.

CardioSoft Initial Screen


Assembly and Setup

NOTE:

Depending on the version of software and user preferences, the initial screen will vary.

7. On the upper right side of the screen, click the **System Configuration** button.

Depending on the version of software the button may look like one of the following:

The **System Configuration** window opens.

8. Click the **Devices** tab.
9. The COM port number is displayed in the **Port** field, as shown in the following screen.

10. Ensure the **Treadmill Model** is set to **T2100**.
11. Ensure the COM port number is set to the port number as identified in step 5.d, above.

If necessary, change the number to match.

Assembly and Setup


12. Click **OK** to close the **System Configuration** window and return to the CardioSoft or CS initial screen.

If needed, the CardioSoft/CS operating instructions should be referenced for controlling the treadmill with the CASE application.

MAC 5500 ST connection and setup:


NOTE:

Prior to starting this section, confirm the MAC 5500 device contains the ST (Stress Test System) option. The device can be identified by the two additional rows of keys on the keyboard, as shown (1) in the following illustration.

MAC 5500 ST (Keyboard)

1. Turn the power **OFF** at the treadmill.
2. Connect the 9-pin D-Sub end of the communication cable (GE part number 2097829-108) to the RS232 COM port located on the lower front panel of the treadmill.
3. Connect the 8-pin mini DIN of the communication cable (GE part number 2097829-108) to the COM1 port on the MAC 5500 ST unit.

MAC 5500 Rear View

4. Make sure the treadmill is connected to AC power.
5. Turned the power **ON** at the treadmill.
6. If necessary, turn on the MAC 5500 ST device with the power button on the keyboard.


Assembly and Setup

Refer to the MAC 5500 operator's manual for instructions on conducting an exercise stress test.

MAC 2000 ST connection and setup:


NOTE:

Prior to starting this section, confirm that the MAC 2000 device contains the ST (Stress Test System) option. The device can be identified by the two additional rows of keys on the keyboard, as shown (1) in the following illustration.

MAC 2000 ST (Keyboard)

1. Turn the power **OFF** at the treadmill.
2. Connect one end of the RS232 serial cable (GE part number 2097829-110) to the RS232 COM port located on the lower front panel of the treadmill.
3. Connect the other end of the RS232 serial cable (GE part number 2097829-110) to the COM B port on the MAC 2000 ST unit.

MAC 2000 Rear View

4. Make sure the treadmill is connected to AC power.
5. Turn the power **ON** at the treadmill.
6. If necessary, turn on MAC 2000 ST with the power button on the keyboard.

Assembly and Setup

Refer to the MAC 2000 operator's manual for instructions on conducting an exercise stress test.

4

Preventive Maintenance

Regular cleaning and maintenance are essential to keep your T2100-ST1 and T2100-ST2 treadmills operating at their best for many years. We recommend that you record all maintenance and service in a log (as shown in Appendix A).

CAUTION

Before cleaning the T2100-ST1 and T2100-ST2 treadmills, turn the main power switch to OFF, and disconnect the treadmill from its power outlet. Never use wet cleaning materials near a power source.

To preserve the condition of your warranty, make sure that all repair procedures (other than normal maintenance) are performed by an authorized and qualified service provider. Contact your local GE customer support representative.

Use only authorized replacement parts. Using other parts may void your warranty and may cause your T2100-ST1 or T2100-ST2 treadmill to malfunction.

NOTE:

Replacement Parts need to be ordered through the manufacturer of all operating equipment to ensure BASIC SAFETY and ESSENTIAL PERFORMANCE with regard to ELECTROMAGNETIC DISTURBANCES for the EXPECTED SERVICE LIFE.

Reprocessing Instructions

To reprocess the treadmill between each patient use, follow the reprocessing instructions in Appendix C.

Preventive Maintenance

Daily Maintenance

- Ensure that the treadmill is functioning properly.
- Visually inspect the treadmill and walking belt for damage and wear.
- Press the Emergency STOP Button (ESB) and the treadmill will coast to a stop. To release the emergency stop switch, turn the push button 1/4-turn in counterclockwise direction. The treadmill will return to 0.0% elevation.
- Pull the Stop Tether Switch (STS) Switch and the treadmill will have a controlled stop. To re-attach the stop tether, attach to the original position on the switch. The treadmill will return to 0.0% elevation.

Weekly Maintenance

- Vacuum around and under the treadmill. Clean all exposed surfaces with a vacuum cleaner. Avoid moving the treadmill from its original position as moving it will compromise the original belt tracking setting.
- Check running belt tension.
- Observe running belt tracking and correct as required.

Monthly Maintenance

Inspect and clean the running belt.

Semiannual Maintenance

- Evaluate the condition of the deck and running belt.
- Adjust the running belt to assure proper alignment.
- Check the running belt adjustment.
- Check the drive belt tension adjustment.
- Clean and lubricate the treadmill elevation screw.
- Clean the interior of the motor electrical enclosure as needed.

Preventive Maintenance

Belt Cleaning and Inspection

1. Turn the treadmill main power switch ON.
2. Start treadmill at 0.5 mph (0.8 km/h).
3. With a damp small towel wipe excessive dirt from running belt keeping the towel in the center of the length of the treadmill. Avoid getting the towel near the rear roller.
4. When belt is clean, stop the treadmill.
5. Inspect the running belt for tears or nicks. If damaged, replace the belt.
6. Perform the Running Belt Tracking Adjustment and Belt Tension Adjustment.

Running Belt Tracking Adjustment

This procedure requires the following tool:

- 1/4" Allen wrench

NOTE:

Because this adjustment is not covered under your warranty, it is important that you review these instructions thoroughly before proceeding.

The patented MasterTrack® Belt Tracking System significantly reduces the need to adjust the belt on your T2100-ST1 or T2100-ST2 treadmill. However, when you operate your T2100-ST1 or T2100-ST2 treadmill for the first time, you may need to adjust the tracking of the belt to conform to your floor. You may also need to adjust the tracking if you move the machine to another location.

Your running belt should remain centered, although a slight amount of movement to the left or right is normal during use. Do not allow the running belt to travel all the way to either side.

To adjust the belt tracking, do the following:

1. Turn the treadmill's power switch to ON.
2. Increase the speed to 3 mph (4.8 km/h)
3. Observe the left side of the running belt as it travels over the rear roller. If the belt runs to the right side of the roller, turn the right bolt 1/8 turn clockwise, and turn the left bolt 1/8 turn counterclockwise.

NOTE:

When tightening one side of the belt, always loosen the opposite side one-half as much. This procedure provides finer control, with a smaller impact on the belt tension.


Check the belt after 2 minutes, with the treadmill running at approximately 7 mph (11.3 km/h). If the belt does not correct

Preventive Maintenance

itself, continue with slight turns until the belt is in the center of the rear roller. If the belt runs toward the left side of the roller, reverse the adjustments.

NOTE:

Uneven floors accelerate belt misalignment. This situation may require more frequent adjustment to prevent belt damage.

Item	Description
1	Adjustment Bolts

NOTE:

When operating the treadmill at a High-Speed Application may cause hesitation/slippage of the running belt with each foot plant. This could be a sign of the backing of the running belt breaking down causing a premature failure. For troubleshooting refer to "Running Belt High Speed Application" on page 112.

Running Belt Tension Adjustment

The running belt may stretch and loosen with regular use. This looseness is noticeable when the belt tends to hesitate or stick. Adjust the tension on the belt by following the procedure and referring to the illustration above.

1. Turn the treadmill's power switch to ON.
2. Start the treadmill and increase to 1 mph (1.6 km/h).
3. Start walking on the treadmill, grabbing side handrail and applying pressure with your foot to create resistance on the running belt.
4. If the running belt hesitates or slips on the front drive roller, tighten both tension bolts $\frac{1}{2}$ turn (clockwise).
5. Repeat steps 2 through 4 until the running belt stops slipping.

NOTE:

When the running belt is too tight, the edge of belt will curl, causing premature running belt failure.

Preventive Maintenance

Drive Belt Tension Adjustment

This procedure requires the following tool:

- 1/8" Allen Wrench
- 3/4" Socket or Box Wrench
- 3/4" Wrench
- Tape measure

NOTE:

The drive belt may stretch and loosen with regular use. This looseness may result in a flapping noise under the hood.

If there are indications that the drive belt has stretched and become loose, contact your in-house biomedical technician, or contact GE Medical Systems *Information Technologies, Inc.*, and open a customer service call.


1. Remove the (5) #10-32 screw located on the bottom hood with 1/8" Allen wrench.
2. With the 3/4" socket and wrench loosen the TENSION ADJUSTMENT back nut.
3. Press down on drive belt between motor and front roller with approximately 5lbs of force to achieve 1/4" to 3/8".
4. If adjustment is need using a 3/4" wrench turn the TENSION BOLT clockwise to tighten. To loosen, turn the wrench counterclockwise.
5. Make small adjustment until the drive belt deflects approximately 1/4" to 3/8", tighten the TENSION ADJUSTMENT back nut.

NOTE:

Failure to lock the TENSION ADJUSTMENT back nut will allow the drive belt tension to become loose.

Preventive Maintenance

6. When the treadmill is properly adjusted, reinstall the hood with (5) #8-32 screws.

Item	Description
1	1/4 to 3/8 Deflection
2	5 lbs. down force
3	Tension adjustment back nut

NOTE:

When drive belt is over-tensioned, the belt tension will cause motor noise. This could result in premature motor life failure.

Exterior Care

The powder-coat finish on your T2100-ST1 or T2100-ST2 treadmill is an extremely durable finish and requires minimal care. Do not allow perspiration to build up on your treadmill. Wipe the unit daily.

- Use a moist cloth to wipe the surface clean; do not allow liquids to enter the system. All cleaning agents and disinfectants used in hospitals and containing up to 70% alcohol are suitable. If liquids have entered the system, notify service to have the system inspected for damage before it is used again.
- DO NOT use disinfection with a phenol base or peroxide compound to disinfect the external surface.

Preventive Maintenance

Elevation Screw Lubrication

The Elevation Screw must be cleaned and lubricated every 6 months to maintain proper operation of the treadmill. Failure to perform this maintenance function will result in premature wear and ultimate failure of the lift mechanism.

This procedure requires the following tools:

- Service Grease, 6 oz tube (GE part number 2097829-072)
 - Clean, lint-free cloth
 - Small paint brush
1. Raise the treadmill to its maximum elevation.
 2. Turn the main power switch to the OFF position and unplug the treadmill from its outlet.
 3. Using a lint-free cloth, remove the old lubricant and accumulated dust from the elevation screw.
 4. Use a small brush to reapply a thin coat of grease to the threads of the elevation screw. Do not use too much grease (less than 1 ounce required)—the excess could squeeze onto the floor and create a slip-and-fall hazard.
 5. Return the unit to service.

Running Deck Maintenance

The treadmill's running deck is maintenance-free and offers two running deck surfaces for double the life of ordinary treadmills.

NOTE:

Do not use any silicone or lubricating sprays during preventative maintenance on your treadmill deck. Using such sprays may void the warranty. Using such sprays can bring about surface changes that may result in hesitation or excessive belt slip.

Prolonged use in high-speed running may cause hesitation/slippage on each foot plant. Inspect the running deck for factory lubrication on running surface. For troubleshooting refer to "Running Belt High Speed Application" on page 112.

If the running deck surface becomes grooved due to wear, it can be renewed by flipping the deck to the opposite side. Notify your in-house biomedical technician, or contact GE Medical Systems *Information Technologies, Inc.*, and open a customer service call.

5

Theory of Operation


The T2100-ST1 / T2100-ST2 Universal Power Supply (UPS) treadmill control provides a means of controlling belt speed and elevation for the treadmill using either an RS-232 serial port or USB port. The communication protocol is based on the industry standard. In addition to the physical interfaces with the controlling system (PC), the motor controller and elevation motor, the UPS also includes an interface to an emergency stop switch and an SD card.

Scope

This chapter applies only to the electronic control for the T2100-ST1 and T2100-ST2 treadmills. Specifications for the mechanical treadmill, elevation motor assembly, belt motor, and motor controller are not covered.


T2100-ST Series Block Diagram

4th Edition Smart Power Supply FGLF0495-3 With Quantum/TruTech


Theory of Operation

4th Edition Smart Power Supply FGLF0495-3

3.1 Edition Smart Power Supply FGLF0495-1

Smart Power Supply (PCB) Overview

When power is applied to the unit, the UPS reads the contents of the non-volatile memory and verifies the calibration contents with a checksum. If the contents are not valid, the calibration constants are set to non-operational values and the unit is put into a state where operation is suspended pending recalibration.

If the contents are valid, the unit enters a startup mode where the system checks to see if it should enter self-test mode (described in the “

Theory of Operation

Self-Test Mode" on page 84) and begins to find the zero elevation point by turning on the elevation motor in the down direction. If motion is not detected for two seconds, the system assumes that the lower elevation limit has been reached. During the "elevation zeroing" time, if a command to start the belt is received from the PC, it is assumed that the PC is not aware of the UPS power cycle and the unit will enter a non-responsive state (waiting for either a communication timeout or a stop command). If the elevation limit is found without receiving a start belt command, the UPS enters the operational mode with a set speed and set elevation of zero.

In operational mode, the unit continuously controls the belt speed output, the belt-enable output, and the outputs to raise or lower the elevation based on the current set points for speed and elevation. Through the use of either communication port, the unit accepts commands to modify the belt state, the belt speed, and the deck elevation. The current state of the control may also be queried through the use of this serial port. As commands are received to change the set point speed or set point grade, they are checked for format and range. Pending approval, the new set points are adopted, and adjustments are made to the state of the outputs.

Speed control is done in a closed loop that monitors the belt speed using two different speed references. The first is a Speed Sensor switch located on the motor side of the belt pulley, and the second is a speed reference supplied by the motor controller. A speed calibration is required to provide a target speed output for a given set point. After the initial target is reached, the speed inputs are constantly monitored to adjust the speed output and to detect over speed, zero speed, or speed sensor comparison fault conditions. In the case of a fault condition, the belt is stopped, and operation is suspended.

Elevation position control is done using relays to turn on the elevation motor in the up or down direction depending on the current set point. A micro-switch input is used to determine the position of the elevation motor by counting revolutions of the elevation motor's internal gearing. An elevation calibration is required to provide an accurate number of counts from zero to full elevation. No closed loop control is done for positioning elevation.

Power may be removed from the unit by opening the normally closed switch connected to the emergency stop input, or by removing AC power from the system. In both cases, the outputs will be put into a non-energized state that will correspond to the belt being shut off and the elevation motor being disengaged.

Software Requirements

The following sections describe the software commands that the hosts system uses to control the T2100-ST series' speed and elevation.

Speed Control

The host will send a speed command "SXXX0<CR>" to the processor board. The firmware in the micro controller will translate this command into a data frame, enable the drive if it is not already enabled at this moment, and send the frame through its drive serial interface. The drive will interpret the data frame and will answer with data ACK or NACK, depending on the transmission condition. The speed range is 0.0 to 15.0 mph (0.1 to 24.14 km/h).

Two outputs are required for speed control. The belt enable output must be energized to put the belt into a run mode. The speed control output is a 4kHz 1024 count pulse width modulated (PWM) output supplied directly to the motor controller. Based on calibration, the speed is controlled by outputting a PWM speed reference to the motor controller and enabling the belt run output. The speed inputs are averaged over a two-second period to provide measured speed accuracy of 0.01mph or less.

If the speed of either input is measured to be more than 0.5mph over the set point, or the difference between speed inputs is measured to be more than 0.5mph, the belt will shut off automatically. If the speed measured is more than 3mph over the set point, the belt enable signal will be turned off immediately. If the motor controller activates a fault signal to the UPS, the system will shut down immediately.

Theory of Operation

If the measured speed is within range, it is used to make minor adjustments to the speed output when it is determined that the actual speed and measured speed differ by more than 0.08 mph. These small adjustments (1:1024) are only allowed to be at 4 second intervals to provide for the settle time of the motor and heavily averaged speed input.

Elevation Control

The elevation control process works as follows: The host will send an elevation command "GXXX0<CR>" to the processor board. The micro controller will interpret this data in order to enable the elevation subsystem by asserting EN_ELEV signal on J11 and, based on the actual elevation status, assert/de-assert the UP/DOWN signal on J11 until the requested elevation is achieved. The elevation range is from 0 to 25% grade.

The elevation position is constantly monitored by counting and recording changes in the state of the elevation micro switch sensor. Each debounced state is counted as positive if the elevation motor is engaged in the up direction and negative if it is engaged in the down direction. Using this count and the calibration numbers, the elevation position is controlled by turning on the motor in the desired direction and counting until the target number of counts has been reached. At that time, the elevation motor is disengaged until another set point position is requested.

ESTOP

The emergency stop switch will remove power from the treadmill control board to guarantee a stop condition. When this occurs and the emergency stop button is subsequently released, the control will be in the same state as when the unit is powered up using the main power switch. Therefore, it is necessary in either power up state to guard against inadvertent restarts if the host doesn't recognize the situation and is continuously sending out belt start commands. The method used to guard against this situation is to reject any start command received within 3 seconds of powering up and finding the zero-elevation position.

NOTE:

Allow the treadmill to complete its "Power Up Sequence". If you press the "Start tmill" button on the Host system before the treadmill reaches zero elevation, the treadmill will reject the start command and enters into an error mode. In this error mode, the treadmill will be non-responsive to all start commands. You may observe the elevation to change but belt will not start. You need

to perform a main power cycle on the treadmill to come out of the error mode. During the "Power Up Sequence" there could be a maximum delay of 45 seconds. Pressing the ESB or Stop Tether activation will also serve the purpose of entering the "Power Up Sequence". Alternatively, you may recover from a reject start command error mode by pressing the "STOP tmill" followed by "Start tmill" button.

Communication Hardware (RS-232 Configuration)

- DB9 Connector 9 pin PC/AT Style

Pin	Description
2	Transmit Data
3	Receive Data
5	Signal Ground

- RS232 (+/- 10V)
- 9600 Baud
- No Parity
- 8 Data Bits
- 1 Stop Bit
- Full Duplex

Communication Hardware (USB Configuration)

- PC USB Host Port
- UPS USB Peripheral Port Mini-B Connector
- FTDI FT232RL Virtual Communication Port Driver required
 - 9600 Baud
 - No Parity
 - 8 Data Bits
 - 1 Stop Bit
 - Full Duplex

Theory of Operation

Self-Test Mode

On power up, the UPS sends out a special test character 'T' and checks to see if it is returned. This is a check for a shorting jumper across the transmit and receive lines of the RS-232 serial port. Note that it is not possible to enter self-test using the USB port. If the RS-232 lines are determined to be shorted, the UPS enters the self-test mode. During the self-test mode, 14 stages of speed and elevation set points are used to exercise the inputs and outputs. Upon completion of each new set point, another character is sent out and the process is repeated to advance to the next stage. In this way, the shorting jumper can be removed at any point to stop the staging process. The following table shows the sequence of self-test stages:

Stage	Incline	Speed (mph)
1	5%	0.0
2	10%	0.0
3	15%	0.0
4	20%	0.0
5	20%	1.0
6	20%	2.5
7	20%	5.0
8	20%	7.5
9	20%	10.0
10	15%	7.5
11	10%	5.0
12	5%	2.5
13	0%	1.0
14	0%	0.0

Calibration

Both a speed and an elevation calibration method are provided in the UPS software. The elevation calibration requires manually adjusting the deck to find the sensor counts for 7 fixed elevation positions and is not recommended for field adjustment without special treadmill interface software. During speed calibration, the speed reference output is swept across its range and speed measurements are made at 15 points. This function is only available through the use of special treadmill interface software.

DIP Switch Settings

A dip switch is located on the UPS at location SW2. An arrow indicates the "ON" position, and each switch is numbered 1-4. Currently, only two of the switches are read by the UPS, and both are only checked on power up so changing the switches while the system is running will not have an effect until power is cycled. Switches 1, 2, and 4 are in the "OFF" position, and switch 3 is in the "ON" position. This configuration is compatible with the GE Protocol (Firmware Version 1.57 and later FGLF0495-1 Smart Power Supply).

Additional dip switch on UPS (Firmware Version 1.1812, 1.19.109 and later FGLF0495-03 Smart Power Supply. This configuration has 1-6 number switch. Switch 1,2,4,5 and 6 are in the "OFF" position and switch 3 is in the "ON" position. When retrieving system flash log file switch configuration 1,2,3,4 and 6 are in the "OFF" position and switch 5 is in the "ON" position. Refer to "System Flash Log Retrieval" Troubleshooting section under Chapter-6 Troubleshooting.

Electrical Inputs

- A. *Instrument Power*: 230VAC or 120VAC supplied through two 1/4-amp breakers in parallel for control board power. 230VAC or 120VAC supplied through two 3-amp breakers in parallel for elevation motor power.
- B. *Elevation Sensor Input*: A two-wire interface for detecting the presence of a closed micro-switch on the elevation motor.
- C. *Speed Sensor Input*: A three-wire interface for detecting the state of a hall effect sensor on the motor shaft (only two conductors are used).
- D. *Serial Port Input*: An electrically isolated RS-232 serial port from the host controller.
- E. *USB Port Input*: An electrically isolated USB port from the host controller.
- F. *Motor Controller Fault*: Fault signal supplied from motor controller.
- G. *Voltage Selection Switch*: The onboard input voltage selection switch allows the unit to be factory set to 115VAC or 230VAC.
- H. *Emergency Stop Switch*: A normally closed switch input used to switch power to the UPS control must be rated at 24 Volts DC and 3 Amps continuous duty.

Theory of Operation

Electrical Outputs

- A. *Elevation Motor:* AC Common is permanently provided to the elevation motor. In addition, the 120/230VAC hot input for motor elevation power is switched through a 10-amp relay to either the up or down input to the elevation motor.
- B. *Motor Controller Drive Enable Output:* An optically isolated output is provided to switch on the start input circuit for the belt drive. This output is to provide no more than 2mA.
- C. *Speed Output:* A PWM signal for speed is provided to the motor controller. The signal is toggled between the motor controller's signal reference voltages (5Vdc). The maximum current for this output is 2mA and the maximum voltage is 15Vdc.
- D. *Serial Port Output:* An electrically isolated RS-232 level (not to exceed +/- 10VDC) is provided to the host controller.
- E. *USB Port Output:* An electrically isolated USB port to the host controller not to exceed 5V or 100mA.

Electrical Connections

- A. *Input Power Header:* An amp - 641968-1 four-pin MATE & LOCK header is provided for instrument power.
- B. *Elevation Output Header:* An amp 641968-1 three-pin MATE & LOCK headers are provided for output power to the elevation motor.
- C. *Elevation Input Header:* A 3 pin 0.156" pitch header is provided for the input from the elevation sensor.
- D. *Speed Input Header:* A 2-pin 0.156" pitch header is provided for the input from the magnetic reed switch speed sensor.
- E. *Motor Controller Interface:* A 10-pin 0.1" pitch header is provided for interfacing to the motor controller.
- F. *RS-232 Interface:* A DB9 Connector 9-pin PC/AT Style is provided for interfacing with the host controller.
- G. *USB Interface:* A Mini-B right angle USB Connector is provided for interfacing with the host controller.
- H. *Emergency Stop Switch:* Two 0.187" quick-connect blade terminals are provided for an emergency stop switch input.

Physical Requirements and Restrictions

- A. The printed circuit board should be a rigid 1/16" thick, FR4, SMOBC circuit board with 2-ounce copper.
- B. Galvanic isolation between the AC power circuits and the low voltage electronic circuits should be guaranteed.
- C. Outer Dimensions: 6" T x 9" W x 2.5" D (max).
- D. *Mounting position:* The board shall be capable of being mounted with the PCB either parallel or perpendicular to the bottom pan of the treadmill.
- E. *Environment:* 0 - 70 °C ambient temperature.
- F. High level vibrations from adjacent equipment may be present.
- G. After installation into the treadmill, expect the board to experience moderate shock impulses during shipment.

6

Troubleshooting

The T2100-ST Series Calibration

Software is designed to assist in gathering information when troubleshooting the T2100 ST1 and T2100-ST2 treadmill.

For information on how to connect the T2100-ST Series Calibration Software and retrieve error codes, refer to "Calibration" starting on page 199. A description of the "**Error Code Tab**" can be found on page 217.

When you receive an error code from the T2100-ST series treadmill, use the flowcharts in this chapter to help you identify and resolve the error. Descriptions of the error codes are provided in this chapter.

There are two methods for using the flowcharts.

If you know the failure mode causing the issue, scan the flowchart titles for that failure mode and review the corresponding flowchart.

If you are receiving an error code, locate the error code in the "Troubleshooting Guidance Base on Error Code" table, below, and then refer to the referenced flowchart(s).

If you still cannot identify or resolve the error, contact local GE customer support representative.


Troubleshooting Guidance Base on Error Code

Error Code Reference	Possible Failure Mode / Possible Root Cause	Related Troubleshooting Guidance
01, 08	Running belt will not move / elevation reaches target	- Smart Power Supply Error Code 1 Flow Chart 1J "Bad Calibration Error (1)" on page 100. - Smart Power Supply Error Code 8 Flow Chart 1L "Motor Controller Fault Signal (8)" on page 109.
02	Treadmill would not elevate / "No power to actuator"	- Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)" on page 101. - Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)" on page 102
02	Treadmill would not elevate and return to park position / "Not receiving count"	- Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)" on page 101. - Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)" on page 102
02	Treadmill will increase elevation but will not decrease elevation / "Smart Power Supply Relay Activation"	- Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)" on page 101. - Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)" on page 102
02	Treadmill would not elevate under heavy load.	- Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)" on page 101. - Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)" on page 102
03, 05	Treadmill will not increase speed during stress test.	- Smart Power Supply Error Code 3 Flow Chart 1L "Zero Speed Motor Controller Error (3)" on page 103. - Smart Power Supply Error Code 5 Flow Chart 1N "Zero Speed External Sensor Error (5)" on page 105.
04, 06	Treadmill speed will not increase above minimum speed.	- Smart Power Supply Error Code 4 Flow Chart 1M "Over Speed Motor Controller Error (4)" on page 104. - Smart Power Supply Error Code 6 Flow Chart 1O "Over Speed External Sensor Error (6)" on page 106.

Troubleshooting


Error Code Reference	Possible Failure Mode / Possible Root Cause	Related Troubleshooting Guidance
04	Treadmill running belt shutdowns around 2 mph.	Smart Power Supply Error Code 4 Flow Chart 1M "Over Speed Motor Controller Error (4)" on page 104.
06, 07	Treadmill running belt shutdown above 6 mph.	- Smart Power Supply Error Code 6 Flow Chart 1O "Over Speed External Sensor Error (6)" on page 106. - Smart Power Supply Error Code 7 Flow Chart 1P "Speed Compare Error (7)" on page 107.
09	Treadmill will start belt after cycling the Emergency Stop Button.	Smart Power Supply Error Code 9 Flow Chart 1R "Belt Start Reject Error (9)" on page 110.
10	Treadmill running belt shutdown due to uncommanded stop/acceleration	Consult the factory for trouble shooting.

Incoming Power 110-240VAC Flow Chart 1A


Troubleshooting

Incoming Power Inline Filter Flow Chart 1B


Incoming Power Drive PC2303-012-N Flow Chart 1C


Troubleshooting

Smart Power Supply Incoming Power 110-240VAC Flow Chart 1D


Emergency Stop Flow Chart 1E


Troubleshooting

Stop Tether Switch Flow Chart 1F


Communication RS232 Flow Chart 1G


Troubleshooting

Communication USB Flow Chart 1H


Smart Power Supply Error Code Identification Flow Chart 1I


Troubleshooting

Smart Power Supply Error Code 1 Flow Chart 1J "Bad Calibration Error (1)"

Smart Power Supply Error Code 1 Flow Chart 1J


Smart Power Supply FGLF0495-1 Error Code 2 Flow Chart 1K "Elevation Error (2)"


Troubleshooting

Smart Power Supply FGLF0495-03 Error Code 2 Flow Chart 1KK "Elevation Error (2)"


Smart Power Supply Error Code 3 Flow Chart 1L "Zero Speed Motor Controller Error (3)"


Troubleshooting

Smart Power Supply Error Code 4 Flow Chart 1M "Over Speed Motor Controller Error (4)"


Smart Power Supply Error Code 5 Flow Chart 1N "Zero Speed External Sensor Error (5)"

Troubleshooting


Smart Power Supply Error Code 6 Flow Chart 10 "Over Speed External Sensor Error (6)"


```

graph TD
    Start([1P]) --> D1{Is the Red LED D13 flashing 7 times with 4 second pause on the Smart Power Supply?}
    D1 -- Yes --> S1[Turn the treadmill power switch OFF (O).]
    S1 --> D2{Was the harness damaged?}
    D2 -- Yes --> R1[Replace Harness]
    D2 -- No --> D3{Was the harness properly connected?}
    D3 -- No --> C1[Connect harness.]
    D3 -- Yes --> D4{Was the harness damaged?}
    D4 -- Yes --> R2[Replace Motor]
    D4 -- No --> D5{Was the harness properly connected?}
    D5 -- No --> C2[Connect harness.]
    D5 -- Yes --> EndB((B))
    D1 -- No --> D2
    D2 --> S2[Verify the 10 pin harness between Drive and Smart Power Supply is properly connected and not damaged.]
    S2 --> I1[Kollmorgen]
    S2 --> I2[Quantum/TruTech]
    I1 --> S3[Verify the harness between Motor and Drive is properly connected and not damaged.]
    S3 --> I3[Kollmorgen 6 Pin Connection]
    S3 --> I4[Quantum/TruTech 11 Pin Connection]
    I3 --> S4[Turn the treadmill power switch OFF (O) for 10 seconds then turn the power switch ON (-).]
    S4 --> D6{Is the Red LED D13 on the Smart Power Supply on continuously?}
    D6 -- No --> EndA([Go To Page 1I])
    D6 -- Yes --> EndA
    
```

Troubleshooting


Smart Power Supply Error Code 8 Flow Chart 1L "Motor Controller Fault Signal (8)"

Troubleshooting

Smart Power Supply Error Code 9 Flow Chart 1R "Belt Start Reject Error (9)"

NOTE:

T2100-ST using Firmware Version 1.18.12 and previous. During power up sequence the treadmill communication is maintained. Pressing the "Start tmill" button before completing Power Up Sequence, will result in Error Code 9.

T2100-ST using Firmware Version 1.19.109 and later. During power up sequence the treadmill does not establish communication until Power Up Sequence is complete. This prevents the generation of Error Code 9 even if "Start tmill" is pressed during Power Up Sequence.

Kollmorgen Drive PC2303-012-N status LED CODE list

When the drive is faulted, the Status LED blinks out the fault code. The Status LED blinks the proper number of times, pauses with the Status LED off, and then repeats the blink code.

Status Code	Description	Fault Retention
No Blink	No faults, power stage enabled	-
Fast Blink	No faults, power stage disabled	-
LED Off	Ac line power insufficient	-
1 Blink	PWM_CMD Stuck Low	Latched (See Note 1)
2 Blink	Gcoder Feedback Error	Latched (See Note 1)
3 Blink	Not Used	-
4 Blink	Output Over Current	Latched (See Note 1)
5 Blink	Control Supply Under Voltage	Self-Resetting
6 Blink	Drive Over Temperature	Latched (See Note 1)
7 Blink	Bus Over Voltage	Self-Resetting
8 Blink	Output Short Circuit	Latched (See Note 1)

NOTE:

Contact local GE customer support representative for possible Drive PC2303-012-N replacement.

Quantum Drive status CODE list

When the drive is faulted, the Status fault code is reported to the Smart Power Supply Board.

Code	Fault	Description
0	Bus Over-Current	Bus current too high (shorted motor)
1	Bus Undervoltage	Bus voltage too low.
2	Bus Overvoltage	Bus voltage too high.
3	Motor Overload	Excessive load for an extended time.
4	Motor Overspeed	Motor RPM is above the maximum.
5	Motor Stalled	Motor RPM stuck at zero.
6	Excess Brake	Brake (Regen Load) dissipation is too high.
7	HeatSink Overtemp	Heatsink temperature is too high.
8	Encoder Pwr Overload	5V supply to encoder is overloaded (less than 4.5V)
9	Initial Speed Cmd High	PWM speed command input duty is too high when enable goes active.

Troubleshooting

10	Speed Command Err	PWM speed command frequency is too high.
11	Invalid Hall state	Hall sensor malfunction

NOTE:

To retrieve Status fault code see "System Flash Log Retrieval" for Manufacturer review.

Contact local GE customer support representative for possible drive replacement.

Running Belt High Speed Application

When operation the treadmill at a High-Speed Application. The running belt is displaying hesitation/slippage of the running belt with each foot plant. The following may be detected:

- Inspect the backing of the running belt for a fabric breakdown, replace the running belt and flip running deck.
- If grooves are found on the running deck. Inspect the running belt fabric. If running belt backing exhibits abnormal wear replace running and flip running deck.

NOTE:

The running deck is reversible providing a new running surface. If both sides of the running deck have been used replace running deck.


System Flash Log Retrieval

When retrieving on UPS (Software Version 1.1812, 1.19.109 and later FGLF0495-03 Smart Power Supply system flash log file power the treadmill down. Set switch configuration 1,2,3,4 and 6 are in the "OFF" position and switch 5 is in the "ON" position.

Insert a blank micro-SD card into IC3, then powerup the treadmill. Depending on the number of records to be written, it can take up to 2 minutes to finish. When complete the treadmill begins to initialization by turning the down relay on. Once that is done, power the treadmill down and retrieve the micro-SD card. Then reset the switch

configuration 1,2,4,5 and 6 are in the "OFF" position and switch 3 is in the "ON".

NOTE:

The system flash log does not exist on FGLF0495-1 Smart Power Supply.

The micro-SD card cannot be larger than 32 gb with a formatting of FAT32 configuration to download log file.

Released


7

Field Replaceable Units


This chapter lists the parts that can be replaced onsite, identifies their locations, and provides their orderable part numbers.

Field Replaceable Units

Final Assembly

Final Assembly Circuit Board Connection

Orderable FRUs

Item	Description	GE Orderable Number	Qty
2	Handrail RH	2097829-066	1
3	Handrail LH	2097829-067	1
8	Handrail Center	2097829-068	1
10, 12	Emergency Stop with Guard Assembly	2097829-011	1
10, 11	Pull Tether Assembly	2097829-026	1
19, 20, 21, 28, 29	Hood Kit Assembly	2097829-069	1
30	Comm Cable Clamps	2097829-074	1
22, 23, 24	Test Plug	2097829-070	1
4-7, 9, 13, 16-18, 23, 24	Final Assembly Hardware Kit	2097829-071	1

Field Replaceable Units

Item	Description	GE Orderable Number	Qty
14	POWER CORD 125V C19 NEMA 5-20P (USA)	2097829-101	1
	POWER CORD 250V C19 NEMA 6-15P (USA)	2097829-102	
	POWER CORD 250V 16A C19 CEE 7/7 (EU)	2097829-103	
	POWER CORD 250V 13A C19 BS1363 (UK)	2097829-104	
	POWER CORD 250VAC 15A C19 BS546 3 PIN (India)	2097829-105	
	POWER CORD 250VAC 15A C19 AS/NZS 3112 (ANZ)	2097829-106	
	POWER CORD 250VAC 16A C19 GB 1002 (China)	2097829-107	
	POWER CORD 250VAC 13A TYPE K-DANISH	2097829-112	
	POWER CORD 250VAC 16A TYPE J-SWISS	2097829-113	
	POWER CORD 250VAC 16A TYPE H-ISRAEL	2097829-114	
	POWER CORD 250VAC 15A TYPE L-ITALY	2097829-115	
	POWER CORD 250VAC 15A TYPE N-BRAZIL	2097829-116	
	POWER CORD 250VAC 16A BS546 INDIAN	2097829-117	
	POWER CORD 250VAC 20A NEMA 6-20P	2097829-118	
Not shown	Grease-tube	2097829-072	1


For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
2	Handrail RH	1
3	Handrail LH	1
4	Washer Nylon	2
5	Lock Washer 3/8"	4
6	Screw HSHCS 3/8"-16 x 1" Lg PLTD	6
7	Dome Plug	4
8	Handrail Center	1
9	Screw HSBHCS 1/4"-20 x 3/4" Lg	8

Item	Manufacturer's Description	Qty
10	Harness ESB & Pull Tether	2
11	Pull Tether Assembly or	1
12	Emergency Stop Assembly IDEC	1
13	Riv-Nut 1/4"-20	1
14	POWER CORD 125V C19 NEMA 5-20P (USA) POWER CORD 250V C19 NEMA 6-15P (USA) POWER CORD 250V 16A C19 CEE 7/7 (EU) POWER CORD 250V 13A C19 BS1363 (UK) POWER CORD 250VAC 15A C19 BS546 3 PIN (India) POWER CORD 250VAC 15A C19 AS/NZS 3112 (ANZ) POWER CORD 250VAC 16A C19 GB 1002 (China) POWER CORD 250VAC 13A TYPE K-DANISH POWER CORD 250VAC 16A TYPE J-SWISS POWER CORD 250VAC 16A TYPE H-ISREAL POWER CORD 250VAC 15A TYPE L-ITALY POWER CORD 250VAC 15A TYPE N-BRAZIL POWER CORD 250VAC 16A BS546 INDIAN POWER CORD 250VAC 20A NEMA 6-20P	1
16	Clamp Support	1
17	Lock Washer 1/4"	1
18	Bolt HHCS 1/4"-20 x 3/4" Lg Gr 5	1
19	Hood T2100-ST Series with GE 80mm Logo	1
20	Decal Caution	1
21	Screw HHWHTS #8-32 x 3/8" Lg	5
22	Test Plug	1
23	Velcro Hook Backing 1" Lg.	1
24	Velcro Backing 1" Lg.	1
28	DECAL GE LOGO 80 MM	1
29	STOP DECAL	1
30	Comm. Cable Clamp	2
1, 15, 25, 26, 27	Not orderable as FRU	

Field Replaceable Units

Deck Assembly

Orderable FRUs

Item	Description	GE Orderable Number	Qty
2	Belt Running with V-Guide	2097829-057	1
7	Deck Isolator 1.75" Dia with Stud	2097829-058	1
9	Deck Running Board	2097829-059	1
11	Roller Front	2097829-060	1
15	Rear Roller	2097829-061	1
18	End Cap LH	2097829-062	1
19	End Cap RH	2097829-063	1
21	Extrusion Anti-Skid Rubber 60" LG	2097829-064	1
3-6, 8, 10, 12-14, 16, 17, 20	Deck Assembly Hardware Kit	2097829-065	1


For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
1	Not orderable as FRU	
2	Belt Running with V-Guide	1
3	Bolt HHCS 5/16-18 x .75" Lg	6
4	Lock Washer 5/16" Dia	7
5	Flat Washer 5/16" Dia	9
6	Deck Support Channel Weldment	3
7	Deck Isolator 1.75" Dia with Stud	6
8	Flange Lock Nut HNS 5/16-18	6
9	Deck Running Board	1
10	Screw HSFCHCS 5/16"-18 x 1.5" Lg PLTD	6
11	Roller Front	1
12	Set Screw 5/16"-18 x 1.75" Lg	1
13	Nut FHNZ 5/16-18 Grade 5	1
14	Bolt HSHCS 5/16"-18 x 2" Lg PLTD	1
15	Rear Roller	1
16	Belt Tension Plate	2
17	Screw HSHCS 5/16"-18 x 3" Lg PLTD	2
18	End Cap LH	1
19	End Cap RH	1
20	Screw HSBHCS 5/16"-18 x 1" Lg PLTD	2

Field Replaceable Units

Item	Manufacturer's Description	Qty
21	Extrusion Anti-Skid Rubber 60" LG	2

Frame Assembly

**ACTUATOR INSTALLATION NOTES:**

1. PRIOR TO INSTALLATION OF BOLT, TURN ACTUATOR END CLOCKWISE UNTIL IT STOPS. NEXT TURN COUNTERCLOCKWISE 1 1/2 TURN OUT PLUS ANY ADDITIONAL REVOLUTION TO ALLOW BOLT ALIGNMENT.
2. WHEN TIGHTENING ELEVATION BOLT AND NUT, THE WASHERS MUST MOVE FREELY.

Orderable FRUs

Item	Description	GE Orderable Number	Qty
5	Stem Bumper	2097829-052	1
13	Bumper 1.8" Dia Black	2097829-053	1
21	Belt Drive	2097829-054	1
14-17	Landing Gear Pivot Hardware Kit	2097829-055	1
8-11, 18-20, 22	Frame Assembly Hardware Kit	2097829-056	1


For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
5	Stem Bumper	2
8	Bolt HHCS 1/4-20 x .75" Lg Gr 5	26
9	Flat Washer 1/4" Dia	26
10	Lock Washer 1/4" Dia	26
11	Nut HN 1/4-20 Grade 5	26
13	Bumper 1.8" Dia Black	2
14	Shoulder Bolt HSHSS 5/8" Dia x 1" Lg	2
15	Elevation Delrin Spacer	4
16	Flat Washer 1/2" Dia	2
17	Nylon Lock Nut NLN 1/2"-13	2
18	Bolt HHCS 3/8-16 x 3" Lg Grade 5	1
19	Washer Nylon	4
20	Nylon Lock Nut 3/8-16	1
21	Belt Drive	1

Field Replaceable Units

Item	Manufacturer's Description	Qty
22	Anti-Slip Safety Walk "Cut 1.5" x 4"	2
1, 2, 3, 4, 6, 7, 12	Not orderable as FRU	

Motor Pan Assembly Kollmorgen

Orderable FRUs

Item	Description	GE Orderable Number	Qty
7	Actuator 115V or Actuator 115V MSI (for T2100-ST1 only)	2097829-028 or 2097829-088	1
7	Actuator 220V or Actuator 220V MSI (for T2100-ST2 only)	2097829-029 or 2097829-089	1
12	Fan Assembly DC W/Connector	2097829-030	1
15	Resistor Omite 85/150MJ40R0DE	2097829-031	1
18	Filter CorCom F4635 Or, Filter CorCom 10EHT1	2097829-032 2097829-082	1
19, 20	Drive DC Servo PC2303-012-N with Regen Model RA2302	2097829-033	1
24	Breaker 3 amp	2097829-034	1
25	Breaker 1/4 amp	2097829-035	1
27	Switch TA35-C324-L200C0 ALT Switch MCProducts 2421- 305-119-2000	2097829-036 OR 2097829-123	1
30	Power Cord Inlet C-20 4793.4000	2097829-037	1
4, 5, 21, 22	Motor Pan Isolation Kit	2097829-038	1
2, 3, 6, 8-11, 13, 14, 16, 17, 21, 22, 23, 26, 28, 29, 31,	Motor Pan Hardware Kit	2097829-039 OR 2097829-124	1


For Reference - Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
2	HHCS 1/4-20 x 1.25 Lg Gr 2	3
3	FW 1/4" Dia (above) and Fender Washer 1/4" Dia 1" OD (below)	4
4	Motor Isolation Stem	4
5	Motor Isolation Spacer	4
6	NLN 1/4-20	1
7	Actuator 115V (for T2100-ST1 only)	1
	Actuator 220V (for T2100-ST2 only)	1


Field Replaceable Units

Item	Manufacturer's Description	Qty
8	HHCS 3/8-16 c 5.0" Lg Gr 5	1
9	FW 3/8" Dia	3
10	Washer Nylon	3
11	NLN 3/8-16	1
12	Fan Assembly DC W/Connector	1
13	RHMS #6-32 x .50" Lg	4
14	K-Lock Nut #6-32	4
15	Resistor Omite 85/150MJ40R0DE	1
16	HNS #8-32 Flange	7
17	HHWHTS #8-32 X .50" Lg	7
18	Filter CorCom F4635 or 10EHT1 NOTE: F4635 is being replaced by backward compatible 10EHT1 part.	1
19, 20	Drive DC Servo PC2303-012-N with Regen Model RA2302	1
21	RHMS #10-32 X .375 LG	6
22	NYLON ISOLATOR #10-32	5
23	Plate Breaker Trackmaster	1
24	Breaker 3 amp	2
25	Breaker 1/4 amp	2
26	HHWHTS #8-32 x .375" Lg	2
27	Switch TA35-C324-L200C0	1
28	RHMS #6-32 x .50" Lg	2
29	K-Lock Nut #6-32	2
30	POWER CORD INLET C-20 4793.4000	1
31	PLATE SWITCH	1
1.	Not orderable as FRU	

Motor Pan Assembly Kollmorgen Wiring

Field Replaceable Units

Orderable FRUs

Item	Description	GE Orderable Number	Qty
3	Harness 110V DC Treadmill (for T2100-ST1 only)	2097829-040	1
	Harness 220V DC Treadmill (for T2100-ST2 only)	2097829-041	1
4	Harness Resistor	2097829-042	1

Item	Description	GE Orderable Number	Qty
5	Speed Sensor Hall Effect (2 wire) or, Speed Sensor Hall Effect (3 wire)	2097829-043	1
		2097829-075	
	Harness Speed Sensor Hall Effect (2 wire) Or, Harness Speed Sensor Hall Effect (3 wire)	2097829-044 2097829-076	1
6	Harness Elevation FGLF0495-1 Or, Harness Elevation FGLF0495-3	2097829-045 2097829-083	1
8	Harness Motor Controller	2097829-046	1
10	Harness Main DC Treadmill FGLF0495-1 Or, Harness Main DC Treadmill FGLF0495-3	2097829-084 2097829-047	1
13	Surge Suppressor	2097829-048	1
15	Cable RS232 Panel Mount	2097829-049	1
16	Cable USB Panel to USB B	2097829-050	1
1, 2, 7, 9, 11, 12, 17-23	Motor Pan Assembly Wiring Hardware Kit	2097829-051	1
24, 25, 26	Relay 110V (For T2100-ST1 only) includes Velcro hook/loop	2097829-085	1
24, 25, 26	Relay 220V (For T2100-ST2 only) includes Velcro hook/loop	2097829-086	1
25	Velcro Hook x 1" wide (included in the relay kit)	NA	
26	Velcro Loop x 1" wide (included in the relay kit)	NA	

NOTE:

When replacing old elevation or main wire harness with the new ones, replace both item 6 and 10 with the correct voltage relay.


For Reference – Manufacturer's Description for individual items

Item	Manufacturer's Description	Qty
1	Screw HHWTS #8-32 x .375 Lg	4
2	Clamp Plastic 1/4" ID	4

Field Replaceable Units

Item	Manufacturer's Description	Qty
3	Harness 110V DC Treadmill (for T2100-ST1 only)	1
	Harness 220V DC Treadmill (for T2100-ST2 only)	1
4	Harness Resistor	1
5	Speed Sensor Hall Effect	1
	Harness Speed Sensor Hall Effect	1
6	Harness Elevation	1
7	HH Screw #10-32 x .375 Lg GRD	2
8	Harness Motor Controller	1
9	Harness Ground	2
10	Harness Main DC Treadmill	1
11	Harness Jumper	2
12	LW #10 EXT Tooth	1
13	Surge Suppressor	1
15	Cable RS232 Panel Mount	1
16	Cable USB Panel to USB B	1
17	Plate RS232 & USB	1
18	HHCS #8-32 x .50 Lg	2
19	HNS #8-32 Flange Nut	2
20	Decal Caution	2
21	Zip Tie 8" Lg	1
22	Zip Tie 4" Lg	1
23	RMHS #4-40 x .375 Lg	2
14	Not orderable as FRU	
24	Relay 110V (for T2100-ST1 only)	1
	Relay 220V (for T2100-ST2 only)	1
25	Velcro Hook 1" x 2"	1
26	Velcro Loop 1" x 2"	1

Motor Pan Assembly Quantum/TruTech

Field Replaceable Units

Orderable FRUs

Item	Description	GE Orderable Number	Qty
7	Actuator 115V or Actuator 115V MSI (for T2100-ST1 only)	2097829-028 or 2097829-088	1
7	Actuator 220V or Actuator 220V MSI (for T2100-ST2 only)	2097829-029 or 2097829-089	1
12	Fan Assembly DC W/Connector	2097829-030	1
15	Resistor Omite 85/150MJ40R0DE	2097829-031	1
18	Filter CorCom F4635 Or, Filter CorCom 10EHT1	2097829-032 2097829-082	1
	PARTS KIT-MOTOR ASSM TRUTECH	2097829-127	1
	PARTS KIT-DRIVE ASSM QUANTUM	2097829-125	1
	BOARD FG0495-3 UNIVERSAL POWER SUPPLY FW V1.19.109	2097829-126	1
24	Breaker 3 amp	2097829-034	1
25	Breaker 1/4 amp	2097829-035	1
27	Switch TA35-C324-L200C0	2097829-036	1
30	Power Cord Inlet C-20 4793.4000	2097829-037	1
4, 5, 21, 22	Motor Pan Isolation Kit	2097829-038	1
2, 3, 6, 8-11, 13, 14, 16, 17, 21, 22, 23, 26, 28, 29, 31,	Motor Pan Hardware Kit	2097829-039	1


For Reference - Manufacture's Description for individual items


Item	Manufacturer's Description	Qty
2	HHCS 1/4-20 x 1.25 Lg Gr 2	3
3	FW 1/4" Dia (above) and Fender Washer 1/4" Dia 1" OD (below)	4
4	Motor Isolation Stem	4
5	Motor Isolation Spacer	4
6	NLN 1/4-20	1

Item	Manufacturer's Description	Qty
7	Actuator 115V (for T2100-ST1 only)	1
	Actuator 220V (for T2100-ST2 only)	1
8	HHCS 3/8-16 c 5.0" Lg Gr 5	1
9	FW 3/8" Dia	3
10	Washer Nylon	3
11	NLN 3/8-16	1
12	Fan Assembly DC W/Connector	1
13	RHMS #6-32 x .50" Lg	4
14	K-Lock Nut #6-32	4
16	HNS #8-32 Flange	5
17	HHWHTS #8-32 X .50" Lg	4
18	Filter CorCom F4635 or 10EHT1 NOTE: F4635 is being replaced by backward compatible 10EHT1 part.	1
23	Plate Breaker Trackmaster	1
24	Breaker 3 amp	2
25	Breaker 1/4 amp	2
26	HHWHTS #8-32 x .375" Lg	2
27	Switch TA35-C324-L200C0	1
28	RHMS #6-32 x .50" Lg	2
29	K-Lock Nut #6-32	2
30	POWER CORD INLET C-20 4793.4000	1
31	PLATE SWITCH	1
1.	Not orderable as FRU	

Field Replaceable Units

Motor Pan Assembly Quantum/TruTech Wiring

Orderable FRUs

Item	Description	GE Orderable Number	Qty
3	Harness 110V DC Treadmill (for T2100-ST1 only)	2097829-040	1
	Harness 220V DC Treadmill (for T2100-ST2 only)	2097829-041	1

Field Replaceable Units

Item	Description	GE Orderable Number	Qty
4	Harness Resistor	2097829-042	1
5	Speed Sensor Hall Effect (3 wire)	2097829-075	1
	Harness Speed Sensor Hall Effect (3 wire)	2097829-076	1
6	Harness Elevation FGLF0495-1 Or, Harness Elevation FGLF0495-3	2097829-045 2097829-083	1
8	Harness Motor Controller	2097829-046	1
10	Harness Main DC Treadmill FGLF0495-3	2097829-047	1
13	Surge Suppressor	2097829-048	1
15	Cable RS232 Panel Mount	2097829-049	1
16	Cable USB Panel to USB B	2097829-050	1
1, 2, 7, 9, 11, 12, 17-23	Motor Pan Assembly Wiring Hardware Kit	2097829-051	1
24, 25, 26	Relay 110V (For T2100-ST1 only) includes Velcro hook/loop	2097829-085	1
24, 25, 26	Relay 220V (For T2100-ST2 only) includes Velcro hook/loop	2097829-086	1
25	Velcro Hook x 1" wide (included in the relay kit)	NA	
26	Velcro Loop x 1" wide (included in the relay kit)	NA	

NOTE:

When replacing old elevation or main wire harness with the new ones, replace both item 6 and 10 with the correct voltage relay.


For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
1	Screw HHWTS #8-32 x .375 Lg	4
2	Clamp Plastic 1/4" ID	4
3	Harness 110V DC Treadmill (for T2100-ST1 only)	1
	Harness 220V DC Treadmill (for T2100-ST2 only)	1
4	Harness Resistor	1
5	Speed Sensor Hall Effect	1
	Harness Speed Sensor Hall Effect	1

Item	Manufacturer's Description	Qty
6	Harness Elevation	1
7	HH Screw #10-32 x .375 Lg GRD	2
8	Harness Motor Controller	1
9	Harness Ground	2
10	Harness Main DC Treadmill	1
11	Harness Jumper	2
12	LW #10 EXT Tooth	1
13	Surge Suppressor	1
15	Cable RS232 Panel Mount	1
16	Cable USB Panel to USB B	1
17	Plate RS232 & USB	1
18	HHCS #8-32 x .50 Lg	2
19	HNS #8-32 Flange Nut	2
20	Decal Caution	2
21	Zip Tie 8" Lg	1
22	Zip Tie 4" Lg	1
23	RMHS #4-40 x .375 Lg	2
14	Not orderable as FRU	
24	Relay 110V (for T2100-ST1 only)	1
	Relay 220V (for T2100-ST2 only)	1
25	Velcro Hook 1" x 2"	1
26	Velcro Loop 1" x 2"	1

Field Replaceable Units

Motor Assembly Kollmorgen

NOTES:

1. WHEN INSTALLING THE FLYWHEEL & DRIVE SPROCKET, TIGHTEN SET SCREWS ON KEYWAY FIRST, AND THEN TIGHTEN THE SECOND SET SCREW TO MOTOR SHAFT.
2. WHEN TIGHTENING ITEMS 5 & 6, CONTINUE TIGHTENING UNTIL MOTOR CANNOT MOVE SIDE TO SIDE.

Orderable FRUs

Item	Description	GE Orderable Number	Qty
1-7	DC Motor Assembly	2097829-020	1
5, 6, 7, 9	Adjustment Hardware Kit	2097829-021	1

For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
1	DC Motor Base	1
2	Motor DC Servo PMCFA-00101-03	1
3	Sprocket Drive 530-5MGT-15	1
4	Flywheel DC W/Magnet .50 dia	1
5	HHCS 7/16-14 x 5" Lg	1
6	NLN 7/16-14	1
7	HHCS 1/2-13 x 3" Lg All Thread	1
9	HN 1/2-13 Gr5	1

Motor Assembly Quantum/TruTech

NOTES:

3. WHEN INSTALLING THE FLYWHEEL & DRIVE SPROCKET, TIGHTEN SET SCREWS ON KEYWAY FIRST, AND THEN TIGHTEN THE SECOND SET SCREW TO MOTOR SHAFT.
4. WHEN TIGHTENING ITEMS 5 & 6, CONTINUE TIGHTENING UNTIL MOTOR CANNOT MOVE SIDE TO SIDE.


Orderable FRUs

Item	Description	GE Orderable Number	Qty
1-9	PARTS KIT-MOTOR ASSM TRUTECH	2097829-127	1
3, 4, 8, 9	Adjustment Hardware Kit	2097829-098	1

For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
1	MOTOR BRACKET TRUTECH	1
2	MOTOR-TRUTECH 53-0481	1
3	CB-7/16-14 X 3.50 GR 2-CLEAR ZINK	1
4	NLN 7/16-14	1
5	MOTOR KEY .188X.188X2.0 LG	1
6	FLYWHEEL TRUTECH	1
7	SPROCKET DRIVE P30-5MGT-15 .625 BORE	1
8	HHCS 1/2-13 3" LONG ALL THREAD	1
9	HN 1/2-13 GR5	1

Drive Board Assembly Quantum/TruTech

Orderable FRUs

Item	Description	GE Orderable Number	Qty
1-7	PARTS KIT-DRIVE ASSM QUANTUM	2097829-125	1


Field Replaceable Units

Item	Description	GE Orderable Number	Qty
3, 4, 5, 6	Drive Board Adjustment Hardware Kit (Quantum)	2097829-099	1

For Reference – Manufacturer's Description for individual items

Item	Manufacturer's Description	Qty
1	DRIVE BRACKET QUANTUM	1
2	DRIVE-QUANTUM 800-0011470R01	1
3	HN #8-32 NYLON	1
4	MOUNT INVERTOR ISOLATION	1
5	HHWHTS #8-32 X .375	1
6	HHCS #8-32 0.50 LG	1
7	DECAL PCB QUANTUM CONTROL	1

Circuit Board Assembly FGLF0495-1

Orderable FRUs


Item	Description	GE Orderable Number	Qty
2	Board FG0495-1 Universal Power Supply	2097829-022	1
3-8	Circuit Board Assembly Hardware Kit for FG0495-1	2097829-023	1

For Reference – Manufacture’s Description for individual items

Item	Manufacturer’s Description	Qty
2	Board FG0495-1 Universal Power Supply	1
3	HN #8-32 Nylon	4
4	Mount Inverter Isolation	4
5	HHWHTS #8-32 x .375	4
6	RHMS #6-32 x 2.5" Lg	1
7	Standoff Round .50 Lg .166 ID	1
8	K-Lock Nut #6-32	1
1, 9	Not orderable as FRU	

Field Replaceable Units

Circuit Board Assembly FGLF0495-3

Orderable FRUs

Item	Description	GE Orderable Number	Qty
2	Board FG0495-3 Universal Power Supply	2097829-080	1
3-6	Circuit Board Assembly Hardware Kit for FG0495-3	2097829-081	1

NOTE:

When replacing FGLF0495-1 universal power supply with FGLF0495-3 new circuit board assembly hardware kit 2097829-081 may be needed.

For Reference – Manufacturer's Description for individual items

Item	Manufacturer's Description	Qty
2	Board FG0495-3 Universal Power Supply	1
3	HN #8-32 Nylon	5
4	Mount Inverter Isolation	5
5	HHWHTS #8-32 x .375	5
1, 6	Not orderable as FRU	

Elevation Assembly

APPLY MOLY BASED GREASE 2X

Orderable FRUs


Item	Description	GE Orderable Number	Qty
3	Wheel Landing Gear	2097829-024	1
2, 4, 5	Elevation Assembly Hardware Kit	2097829-025	1

For Reference – Manufacturer's Description for individual items

Item	Manufacturer's Description	Qty
2	Spanner Wheel	2
3	Wheel Landing Gear	2
4	HHCS 3/8-16 x 3.25 Lg Gr5	2
5	NLN 3/8-16	2
1	Not orderable as FRU	

Field Replaceable Units

Stop Tether Switch Assembly


Orderable FRUs

Item	Description	GE Orderable Number	Qty
1-4	Pull Tether Assembly.	2097829-026	1

For Reference – Manufacture's Description for individual items

Item	Manufacturer's Description	Qty
1	Tether Pull w/Wrist strap	1
2	E-Stop Kill Switch LS500	1
3	Cap 1.50 x 4.00 LPDE BLK Modified	1
4	Harness ESB & Pull Tether (Not shown)	1

Emergency Stop Assembly

Orderable FRUs

Item	Description	GE Orderable Number	Qty
1-6	Emergency Stop with Guard Assembly	2097829-011	1

For Reference – Manufacturer's Description for individual items

Item	Manufacturer's Description	Qty
1	Emergency Stop IDEC	1
2	Cap 1.50 x 4.00 LPDE BLK Modified	1
3	Harness ESB & Pull Tether (Not shown)	1
4	Decal Stop	1
5	RHMS #6-32 x 1/2" Long	1
6	K-Lock Nut #6-32	1

Communication Cables

Orderable FRUs

Item	Description	GE Orderable Number	Qty
1	USB CABLE 5M A TO B (INTERFACE T2100-ST TO CARDIOSOFT/CS)	2044095-001	1
2	CABLE RS232 INTERFACE T2100-ST TO CASE/PC (contains 2097829-009)	2097829-109	1

Field Replaceable Units

Item	Description	GE Orderable Number	Qty
3	SPARE CABLE ANALOG INTERFACE T2100-ST TO MAC5500ST (contains 2097829-008)	2097829-108	1
4	SPARE CABLE RS232 INTERFACE T2100-ST to MAC2000ST (contains 2097829-010)	2097829-110	1

Calibration Software

Orderable FRUs

Item	Description	GE Orderable Number	Qty
1	KIT T2100-ST CALIBRATION SOFTWARE 1_0_1_28	2097829-128	1

Firmware Upgrade

Item	Description	GE Orderable Number	Qty
1	T2100-ST1/ST2 Firmware Upgrade Kit Kolmogorov FV 1.66 & 1.1812	2097829-073	1
	PART KIT-FIRMWARE 1.19.109	2097829-129	1

Optional Accessories

Item	Description	GE Orderable Number	Qty
1	ADJUSTABLE HANDRAIL FOR T2100-ST1 AND T2100-ST2	2097357-003	1
	T2100-ST Tether Cover Guard	2097357-005	1

8

Component/FRU Removal and Replacement

This chapter provides information for removing and replacing field replaceable parts, including required tools, estimated repair duration, and detailed step-by-step instructions.

Established Component Replacement Times

Procedure	Est. Duration
Replacing the Hood	10 min
Replacing the Drive Motor	1hr 30 min
Replacing the Drive	45 min
KM TO QT Motor & Drive Board CONVERSION	2hr 00min
Replacing the Front Roller	40 min
Replacing the Rear Roller	45 min
Replacing the Running Belt	1 hr 30 min
Replacing the Running Deck	30 min
Replacing the Deck Cushion	1 hr 00 min
Replacing the Motor Drive Belt	1 hr 20 min
Replacing/Adjusting the Elevation Actuator	1 hr 30 min
Replacing the Smart Power Supply Relay Board	1 hr 15 min
Replacing the Circuit Breaker	20 min
Relocate Emergency Stop Button (ESB) and Stop Tether Switch	1 hr 00 min
Replacing the Emergency Stop Button (ESB)	30 min
Replacing the Stop Tether Switch (STS)	30 min
Replacing the Right or Left Handrail	30 min
Replacing the Handrail Center	30 min

Component/FRU Removal and Replacement

Procedure	Est. Duration
Removing and reinstalling handrails for moving	45 min
Replacing the Main Power Switch	25 min
Adjusting the Run Belt tracking and tension	25 min
Replacing the Power Cord	10 min
Replacing the CE filter	30 min
Replacing the Ohmite Resistor	30 min
Replacing the Rear Foot	30 min
Replacing the Isolator (6 isolator & belt tracking)	2 hrs 15 min
Replacing the Main Wire Harness FGLF0495-1	30 min
Replacing the Main Wire Harness FGLF0495-3	30 min
Replacing the Elevation Wire Harness FGLF0495-1	20 min
Replacing the Elevation Wire Harness FGLF0494-3	20 min
Replacing the End Cap	15 min
Replacing the Rubber Stripping on Side Channel	15 min
Replacing Both Wheels	10 min
Replacing the RS232 Harness	30 min
Replacing the USB Harness	30 min
Replacing the Relay 110V or 220V	20 min
Recalibrate Speed, includes time to download software	30 min
Recalibrate Speed	20 min
Recalibrate Elevation, included time to download software	40 min
Recalibrate Elevation	30 min

Comprehensive Tool List

Tools:

- Flat head electronics Screwdriver - 3/32" wide
- Flat head Screwdriver – 1/4" wide
- Flat head Screwdriver - 3/8" wide
- T15 Torx screwdriver
- 1/8" Allen Wrench
- 5/32" Allen Wrench
- 1/4" Allen Wrench
- 3/16" Allen Wrench
- 5/16" Allen Wrench
- 3/32" Socket
- 3/16" Open Ended Wrench and Socket

- 1/4" Open Ended Wrench and Socket
- 5/16" Open Ended Wrench and Socket
- 7/16" Open Ended Wrench and Socket
- 1/2" Open Ended Wrench and Socket
- 9/16" (2X) Open Ended Wrench and Socket
- 3/4" Open Ended Wrench and Socket
- Ratchet
- 6" Ratchet Extension
- (2X) Medium length long-nose Pliers
- Wire Cutters or Utility Knife
- Pliers
- Adjustable Wrench
- Multi Meter

Accessories

- Grounded wrist strap
- Service Grease (PN 2097829-072 tube received with treadmill)
- Tape measure
- Computer or laptop
- Calibration Program (T2100-ST Series Calibration Software)
- Stopwatch or watch with a second hand
- Chalk or some other temporary marking device
- Serial or USB Communication Cable


Replacing the Hood

Tools Required

- 1/4" Hex Socket or Flat Blade screwdriver
- Approximate time to complete: 10 min

1. Turn the treadmill's power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the 5 hood-securing screws (2 on the front and 3 on the back side of the treadmill).
5. Remove the hood and set it aside.
6. Reverse the procedure to reinstall hood.

Component/FRU Removal and Replacement

Replacing the Kollmorgen Drive Motor

Tools Required

- 1/4" Socket or Flat Blade screwdriver
- 7/16" Wrench and Socket
- 9/16" Wrench and Socket
- 3/4" Wrench and Socket
- 3/16" Wrench or Socket
- 3/32" Wrench or Socket
- Grounding Wrist Strap

Approximate time to complete: 1. hr and 30 min


1. Elevate the treadmill with the test plug to approximately 15% grade.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood securing screws.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Remove the Motor Feedback 6 pin cable on J1 from the Drive.
8. Remove the Motor Power 4 pin cable on TB2 from the Drive.


9. With the 3/4" wrench and socket, loosen the belt tension and slide the belt off the sprocket, allowing the motor to pivot, gaining access to the motor mount bolts.
10. With the 1/4" hex socket or flat blade screwdriver, remove the (4) 1/4" plastic clamps on the Motor Feedback & Power wires.
11. With the 7/16" wrench and socket, remove the 4 bolts securing the motor mount to the frame.
12. Lift the Motor Assembly out of the treadmill.
13. With the 9/16" wrench and socket, remove the pivot bolt.
14. Remove drive sprocket from the motor shaft by loosening the two 3/32" set screws on its hub, then slide it off.
15. Remove the flywheel from the motor shaft by loosening the two 3/16" set screws on its hub, then slide it off.
16. Slide the flywheel and 1/8" key on the new motor shaft.
Make sure the flywheel does not rub the plastic housing, creating a motor noise.
17. Slide the drive sprocket and 1/8" key on the motor shaft approximately in the same location from old motor.
18. Tighten the setscrews on the hub, first the one with the key, then the other one.
19. Install the new motor to the existing motor base plate.
20. Tighten the pivot bolt until free movement of the motor is removed.
21. Reinstall the motor, making sure the isolator mounts are in place.
When tightening the 4 bolts, tighten equally without over tightening.
22. Adjust the speed sensor to the magnet on the flywheel to approximately 1/4".

23. Reinstall (4) 1/4" plastic clamps on the Motor Feedback & Power wires.
24. Reconnect the Motor Feedback 6 pin cable to J1 on the drive.

Component/FRU Removal and Replacement

25. Reconnect the Motor Power 4 pin cable to TB2 on the drive.
26. Adjust the tension the drive belt. See "Drive Belt Tension Adjustment" on page 73.
27. After installation, run the treadmill and check for excess vibration and noise from the motor area.
28. Calibrate the Belt Speed.
See "Calibrating Speed" on page 206.
29. Replace the hood and (5) screws.

Replacing the Quantum/TruTech Drive Motor Assy

Tools Required

- 1/4" Socket or Flat Blade Screwdriver
- 7/16" wrench and socket
- 5/8" wrench and socket
- 3/4" wrench and socket
- Grounding Wrist Strap


Approximate time to complete: 1. hr and 30 min

1. Elevate the treadmill with the test plug to approximately 15% grade.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood securing screws.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Detach the External Speed Sensor w/ Harness from the Motor Velcro mount. Carefully set Sensor aside.
8. Remove the Motor Power 4-pin cable on J2 from the Drive Board.

9. With the 1/4" hex socket or flat blade screwdriver, remove the (2) 1/4" plastic clamps on the Motor Power 4-pin cable. Remove Zip Tie that secures 4-Pin cable to Motor Pan frame.

Component/FRU Removal and Replacement

10. Remove the Motor Feedback 11-pin cable on J3 from the Drive Board.
11. Detach the Motor Feedback 11-wire cable from the Motor Encoder and remove it from the Motor Pan routing holes.

NOTE:


If the 11-wire cable is not damaged, then skip step 10 and 11.

12. With the 3/4" wrench and socket, remove bottom nut completely to allow the Motor to pivot and gain access to the (4) motor mount bolts.
13. With the 7/16" wrench and socket, remove the 4 bolts securing the motor mount to the frame.
14. Lift the Motor Assembly out of the Motor Pan.
15. Install the new Motor Assembly in the Motor Pan.

NOTE:


Make sure the 4 isolator mounts are in place. When tightening the 4 bolts, tighten equally without over tightening.

16. Attach the External Speed Sensor using Motor Velcro mount. Set gap to magnet on the flywheel to approximately 1/4".

17. Route the Motor Power 4-pin cable and connect to J2 on the Drive Board.
 18. Reinstall the (2) 1/4" plastic clamps on the Motor Power cable and reinstall Zip Tie to secure 4-pin cable to Motor Pan frame.
 19. Route the Motor Feedback 11-pin cable through Motor Pan routing holes.
 20. Reconnect the Motor Feedback 11-pin cable to J3 on the Drive, and to the Motor Encoder.
 21. Install the drive belt and adjust tension. See "Drive Belt Tension Adjustment" on page 73.
 22. After installation, run the treadmill and check for excess vibration and noise from the motor area.
 23. Calibrate the Belt Speed.
- See "Calibrating Speed" on page 206.

24. Replace the hood and (5) screws.

Replacing the Kollmorgen Drive Board Assy

Tools Required

- 1/4" Socket or Flat Blade Screwdriver
- 5/16" Socket with 6" Extension or Nut Drive
- Grounding Wrist Strap

Approximate time to complete: 45 min

1. Elevate the treadmill to approximately 15% grade.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood-securing bolts.

Component/FRU Removal and Replacement

6. Remove the hood and set it aside.

CAUTION:


Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Remove the Motor Feedback 6 pin cable on J1 from the drive.
8. Remove the 10-pin communication cable on J2 from the drive.
9. Remove the Incoming Power 4-pin connector on TB1 from the drive.
10. Remove the Motor Power 4-pin cable on TB2 from the drive.
11. Remove the 3-pin connector from the Regen Circuit board on the front side of the drive.
12. With the 5/16" socket, loosen the front two screws and remove the back screw.

The drive will slide back allowing the drive to be removed.

13. Reverse the procedure to reinstall the drive.
14. Calibrate the Belt Speed.

See "Calibrating Speed" on page 206.

Replacing the Quantum/TruTech Drive Board Assy

Tools Required

- 1/4" Socket or Flat Blade Screwdriver
- 5/16" Socket with 6" Extension or Nut Drive
- 11/32" wrench and socket
- Grounding Wrist Strap


Approximate time to complete: 45 min

1. Elevate the treadmill to approximately 15% grade.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood-securing bolts.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Remove the Incoming Power 4-pin connector on J1 from the Drive Board.

8. Remove the Motor Power 4-pin cable on J2 from the Drive Board.
9. Remove the Motor Feedback 11-pin cable on J3 from the Drive Board.
10. Remove the 10-pin communication cable on J4 from the Drive Board.

Component/FRU Removal and Replacement

11. With the 5/16" socket remove the two front Drive Mount Frame nuts from beneath the Motor Pan and loosen the 1/4" rear screw from above.

NOTE:


With the 11/32" socket secure the two nuts at the bottom of the treadmill and tighten the 1/4" screw at the back with 11/32" nut.

NOTE:

The location of the ground wire under the front of Drive Mount Frame.

12. Lift the front of Drive Board Assembly and slide forward to remove.
13. Reverse the procedure to reinstall the Drive Board Assembly.
14. Calibrate the Belt Speed.

See "Calibrating Speed" on page 206.

Replacing Kollmorgen Motor & Drive Board with Quantum/TruTech Motor & Drive Board Using Conversion Kit

Tools Required

- 1/4" Socket or Flat Blade Screwdriver
 - 7/16" Wrench and Socket
 - 9/16" Wrench and Socket
 - 11/32" Wrench and Socket
 - 5/8" Wrench and Socket
 - 3/4" Wrench and Socket
 - 3/16" Wrench and Socket
 - 3/32" Wrench and Socket
 - Grounding Wrist Strap
- Approximate time to complete: 2 hr and 00 min

1. Elevate the treadmill with the test plug to approximately 15% grade.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood securing screws.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Remove the Motor Feedback 6 pin cable on J1 from the Drive.

8. Remove the Motor Power 4 pin cable on TB2 from the Drive.

Component/FRU Removal and Replacement

9. With the 3/4" wrench and socket, loosen the belt tension and slide the belt off the sprocket, allowing the motor to pivot, gaining access to the motor mount bolts.
10. With the 1/4" hex socket or flat blade screwdriver, remove the (4) 1/4" plastic clamps on the Motor Feedback & Power wires.
11. With the 7/16" wrench and socket, remove the 4 bolts securing the motor mount to the frame.
12. Remove Zip Tie holding motor and feedback cables to frame.
13. Lift the Motor Assembly out of the treadmill.
14. Remove the 10-pin communication cable on J2 from the Drive.
15. Remove the Incoming Power 4-pin connector on TB1 from the Drive.
16. Remove the 3-pin connector from the Regen Circuit board on the front side of the Drive.
17. With the 5/16" socket, loosen the front two screws and remove the back screw.
18. The drive will slide back allowing the Drive to be removed.
19. Remove the 3 nylon standoffs below Drive. Standoffs are not reused and can be discarded.
20. Remove the Ohmite Resistor, two wire harness and nylon standoffs.


NOTE:

This component does not get reused and may be discarded.


21. Remove the old Speed Sensor. Gromet for 2 wire speed sensor needs to be removed before routing feedback cable through.
22. Install the new Drive Board in the Motor Pan.
23. Start the 1/4" rear screw in Motor Pan mounting hole. Slide the rear of the Drive Board Assembly beneath the 1/4" rear screw, and lower front studs through mounting holes in Motor Pan.

NOTE:

Locate the ground wire under the front of Drive Mount Frame.

24. With the 5/16" socket install the two front Drive Mount Frame nuts from beneath the Motor Pan and tighten the 1/4" rear screw from above.
25. Install the new Motor Assembly in the Motor Pan.
Note: Make sure the 4 isolator mounts are in place. When tightening the 4 bolts, tighten bolts on motor to about 2 to 3 threads and do not over tighten.
26. Route the new 3-wire external speed sensor harness, from Smart Power Supply J10 & DC Fan (-) to Motor. Connect to External Speed Sensor and attach the Sensor to Motor using Motor Velcro mount. Set gap to magnet on the flywheel to approximately 1/4".


27. Route the Motor Power 4-pin cable and connect to J2 on the Drive Board.


28. Install the (2) 1/4" plastic clamps on the Motor Power cable and install Zip Tie to secure 4-pin cable to Motor Pan frame.


Component/FRU Removal and Replacement

29. Route the Motor Feedback 11-pin cable through Motor Pan routing holes.


30. Connect the Motor Feedback 11-pin cable to J3 on the Drive, and to the Motor Encoder.

31. Install the drive belt and adjust tension. See "Drive Belt Tension Adjustment" on page 73.

32. Connect the 10-pin communication cable on J4 from the Drive Board. (FIG. A)

(FIG. A)

33. Connect the Motor Feedback 11-pin cable on J3 from the Drive Board. (FIG. A)
34. Connect the Motor Power 4-pin cable on J2 from the Drive Board. (FIG. A)
35. Connect the Incoming Power 4-pin connector on J1 from the Drive Board. (FIG. A)
36. Update the SPSB Firmware to 1.19.109+
37. Calibrate the Belt Speed.
See "Calibrating Speed" on page 206.
38. After calibration, run the treadmill and check for excess vibration and noise from the motor area.
39. Replace the hood and (5) screws.

Component/FRU Removal and Replacement

Replacing the Front Roller

Tools Required

- 1/4" Socket or Flat Blade Screwdriver
- 3/16" Wrench or Socket

Approximate time to complete: 40 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (5) hood-securing bolts.
5. Remove the hood and set it aside.
6. With the 3/16" Allen wrench or socket, remove the button head Allen bolt securing the left and right end caps from the rear of the treadmill.
7. Remove the end caps and set them aside.
8. With the 1/4" Allen wrench or socket, completely remove the left and right tension bolts from the rear roller.

When removing the tension bolts, maintain equal tension on the bolts to avoid stripping the threads on the tension bolt or roller.

9. Push the rear roller forward as far as it will go.
10. With the 1/4" Allen wrench or socket, remove the bolt securing the front roller to the treadmill frame.
11. Push the front roller from the right side of the treadmill to the left side, lift the roller up; slip the timing belt off the front roller, and lift the front roller up and out of the treadmill.

Use caution to avoid pinching the timing belt.

Bolt 5/16-18 x 2" Lg

Lock Washer 5/16"

Flat Washer 5/16"

Front Roller

12. Slide the front roller between the running belt, place the timing belt over the sprocket of the front roller between the running belt, and place the timing belt over the sprocket of the front roller.
13. Replace the front roller.
14. Verify the guide for the belt deck and rollers is to the left during installation.
15. Verify the timing belt goes around the front roller sprocket and is not pinched.
16. Replace the securing bolt and tighten.
17. Align the bolt holes on the rear roller with the access hole of the end caps.
18. Start the tensioning bolts on the rear roller but do not tighten.
19. Replace the left and right end caps and replace the securing screws.
20. Adjust the tension of the running belt.
See "Running Belt Tension Adjustment" on page 72.

Replacing the Rear Roller


Tools Required

- 1/4" Socket or Flat Blade Screwdriver
 - 3/16" Wrench or Socket
- Approximate time to complete: 45 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (5) hood-securing bolts.
5. Remove the hood and set it aside.
6. With the 3/16" Allen wrench or socket, remove the button head Allen bolt securing the left and right end caps from the rear of the treadmill.
7. Remove the end caps and set them aside.
8. With the 1/4" Allen wrench or socket, completely remove the left and right tension bolts from rear roller. When removing the tension bolts maintain equal tension on the bolts to avoid stripping the threads on the tension bolt or roller.

Component/FRU Removal and Replacement

9. Push the rear roller forward as far as it will go.
10. With the 1/4" Allen wrench or socket, remove the bolt securing the front roller to the treadmill frame.
11. Push the front roller from the right side of the treadmill to the left side, lift the roller up; slip the timing belt off the front roller, and lift the front roller up and out of the treadmill.
Use caution to avoid pinching the timing belt.
12. Push the rear roller from the right side of the treadmill to the left side and lift roller up and out to remove it.
13. Replace the rear roller.
14. Slide the front roller between the running belt, place the timing belt over the sprocket of the front roller between the running belt, and place the timing belt over the sprocket of the front roller.
15. Replace the front roller.
16. Verify the guide for the belt deck and rollers is to the left during installation.
17. Verify the timing belt goes around the front roller sprocket and is not pinched.
18. Replace the securing bolt and tighten.

19. Align the bolt holes on the rear roller with the access hole of the end caps.
20. Start the tensioning bolts on the rear roller but do not tighten.
21. Replace the left and right end caps and replace the securing screws.
22. Adjust the tension of the running belt.
See "Running Belt Tension Adjustment" on page 72.

Replacing the Running Belt

Tools Required

- 1/4" Allen Wrench or Socket
- 3/16" Allen Wrench or Socket
- 1/2" Wrench or Socket with Extension
- Jack stand or block capable of supporting the weight of the rear of the treadmill


Approximate time to complete: 1 hr 30 min

1. Elevate the treadmill to 15%.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill plug from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood-securing bolts.
6. Remove the hood and set it aside.
7. Lift the rear of the treadmill and place jack stands or a block to support it.
8. With the 3/16" Allen wrench or socket, remove the bolt securing the left and right end caps from the rear of the treadmill.
9. Remove the end caps and set them aside.
10. With the 1/4" Allen wrench or socket, completely remove the left and right tension bolts from the rear roller.
When removing the tension bolts, maintain equal tension on the bolts to avoid stripping the threads on the tension bolt or roller.
11. Push the rear roller forward as far as it will go.
12. Lift the edge of the running belt up and use the 3/16" Allen wrench or socket to remove the (4) flat head Allen screws that secure the running deck cross braces.
There are (4) screws total, (2) on the left side and (2) on the right side.
13. Push the deck to one side of the treadmill and, from the other side, lift the edge of the deck up and out over the rail and set it aside.
14. With the 1/4" Allen wrench or socket, remove the bolt securing the front roller to the treadmill frame.

Component/FRU Removal and Replacement

15. Push the front roller from the right side of the treadmill to the left side, lift the roller up, slip the timing belt off the front roller and lift the front roller up and out of the treadmill.
Use caution to avoid pinching the belt.
16. Slide the rear roller out of the rear of the treadmill side rails and remove from the treadmill.
17. With the 1/2" wrench or socket with extension, remove the (4) bolts securing the deck cross braces to the treadmill side channels.
18. Using caution to avoid damage to the antistatic tinsel run between the (2) front bolts, remove the cross braces from the treadmill and set the cross braces and the tinsel aside.
19. Verify all pieces of the deck cushions are accounted for and set aside.
20. Remove the old running belt and replace with the new running belt.
21. Verify the guide for the belt deck and rollers is to the left during installation.
22. Replace the (2) deck braces and deck cushions; start the deck braces securing bolts but do not tighten them.
23. Replace the antistatic tinsel using the (2) front deck brace bolts.
24. Replace the rear roller, pushing the rear roller as far forward as possible.
Do not replace or start the tensioning bolts at this time.
25. Verify the guide for the belt deck and rollers is to the left during installation.
26. Slide the front roller between the running belt, place the timing belt over the sprocket of the front roller between the running belt, and place the timing belt over the sprocket of the front roller.
27. Replace the front roller.
28. Verify the guide for the belt deck and rollers is to the left during installation.
29. Verify the timing belt goes around the front roller sprocket and is not pinched.
30. Replace the securing bolt and tighten.
31. Lift the edge of the running belt and slide the deck onto the deck braces.

32. Tighten the cross brace securing bolts.
Use caution to avoid over tightening the bolts.
33. Align the bolt holes on the rear roller with the access hole of the end caps.
34. Start the tensioning bolts on the rear roller but do not tighten.
35. Replace the left and right end caps and replace the securing screws.
36. Adjust the tension of the running belt.
See "Running Belt Tension Adjustment" on page 72.
37. Reinstall the hood.

NOTES:

1. ADJUSTMENT SCREWS MUST BE MATCH TIGHTENED. PROPER TENSION IS ACHIEVED WHEN BELT CAN'T SLIP ON PULLEYS WHEN IN OPERATON.
2. ADJUST DRIVE BELT TENTION PULLING ALL FREE PLAY OUT OF BELT. DO NOT OVER TIGHTEN CAUSING THE BELT TO STRETCH.
3. MOVE REAR PULLEY SIDE TO SIDE IN CHANNELS TO OBTAIN FRONT TO REAR BELT ALIGNMENT AND TRACKING.

Component/FRU Removal and Replacement

Replacing the Running Deck

Tools Required

- 1/4" Allen Wrench or Socket
- 3/16" Allen Wrench or Socket
- 1/2" Wrench or Socket with Extension
- Jack stand or block capable of supporting the weight of the rear of the treadmill

Approximate time to complete: 30 min

1. Elevate the treadmill to 15%.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill plug from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood-securing bolts.
6. Remove the hood and set aside.
7. Lift the rear of the treadmill and place jack stands or a block to support it.
8. With the 3/16" Allen wrench or socket, remove the bolt securing the left and right end caps from the rear of the treadmill.
See the Running Belt Replacement assembly drawing on page 169.
9. Remove the end caps and set aside.
10. With the 1/4" Allen wrench or socket, completely remove the left and right tension bolts from rear roller.
When removing the tension bolts, maintain equal tension on the bolts to avoid stripping the threads on the tension bolt or roller.
11. Push the rear roller forward as far as it will go.
12. Lift the edge of the running belt up and use the 3/16" Allen wrench or socket to remove the six flat head Allen screws that secure the running deck cross braces.
There are (6) screws total: (3) on the left side and (3) on the right side.
13. Push the deck to one side of the treadmill and, from the other side, lift the edge of the deck up and out over the side rail.
The deck is now ready for maintenance, replacement, or to be flipped.

14. Lift the edge of the running belt and slide the deck onto the deck braces.
15. Replace the screws securing the deck to the deck braces and tighten.
16. Verify the guide for the deck is to the left during installation.
17. Tighten the cross brace securing bolts.
Use caution to avoid over tightening the bolts.
18. Align the bolt holes on the rear roller with the access hole of the end caps.
19. Start the tensioning bolts on the rear roller but do not tighten.
20. Replace the left and right end caps and the securing screws.
21. Adjust the tension of the running belt.
See "Running Belt Tension Adjustment" on page 72.
22. Reinstall the hood.

Replacing the Deck Cushion

Tools Required

- 1/4" Allen Wrench or Socket
- 3/16" Allen Wrench or Socket
- 1/2" Wrench or Socket with extension
- Jack stand or block capable of supporting the weight of the rear of the treadmill

Approximate time to complete: 1 hr 00 min

1. Elevate the treadmill to 15%.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill plug from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) hood-securing bolts.
6. Remove the hood and set it aside.
7. Lift the rear of the treadmill and place jack stands or blocks to support it.
8. With the 3/16" Allen wrench or socket, remove the bolt securing the left and right end caps from the rear of the treadmill.
9. Remove the end caps and set them aside.

Component/FRU Removal and Replacement

10. With a 1/4" Allen wrench or socket, completely remove the left and right tension bolts from rear roller.
When removing the tension bolts, maintain equal tension on the bolts to avoid stripping the threads on the tension bolt or roller.
11. Push the rear roller forward as far as it will go.
12. Lift the edge of the running belt up and use the 3/16" Allen wrench or socket to remove the (2) flat head Allen screws that secure the running deck cross brace with the damaged deck cushions.
13. With the 1/2" wrench or socket with extension, remove the (2) bolts securing the deck cross brace to the treadmill side channels.
14. Remove the cross braces from the treadmill.
15. With the 1/2" socket, remove the rubber deck cushion(s) from the brace.
16. Replace the damaged deck cushion and use the 1/2" socket to remove the bolt securing the rubber deck cushion(s) from the brace.
17. Verify the placement of the deck cushion(s).
Refer to the Running Belt Replacement assembly drawing on page 169.
18. Replace the deck brace and deck cushions
Start the deck brace securing bolts, but do not tighten them.
19. Lift the edge of the running belt and slide the deck onto the deck braces.
20. Replace the screws securing the deck to the deck braces and tighten.
21. Verify the guide for the deck is to the left during installation.
22. Tighten the cross brace securing bolts.
Use caution to avoid over tightening the bolts.
23. Align the bolt holes on the rear roller with the access hole of the end caps.
Start the tensioning bolts on the rear roller but do not tighten.
24. Replace the left and right end caps and the securing screws.
25. Adjust the tension of the running belt.
See "Running Belt Tension Adjustment" on page 53.

26. Reinstall the hood.

Replacing the Motor Drive Belt

Tools required

- 1/4" Hex Wrench or Hex Socket
- 3/4" Open-end Wrench & Hex Socket
- 1/4" Hex Nut Driver or Flat Blade screwdriver
- Grounding Wrist Strap
- Tape Measure

Approximate time to complete: 1 hr 20 min


1. Elevate the treadmill to 15%.
2. Turn the treadmill main power switch to the OFF position.
3. Unplug the treadmill plug from the wall receptacle.
4. Allow the treadmill to discharge for 1 minute.
5. Remove the (5) 1/4" hood-securing bolts.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. With the 3/4" wrench and socket, loosen spring tension and the slide belt off the sprocket, allowing the motor to pivot.
8. With the 7/16" socket and wrench, remove the (6) bolts securing the Bracket Motor Right Hand.
9. Loosen the rear belt tension bolt, allowing the rear roller to move forward towards the front roller.
10. Remove the hex socket head cap screw, retaining the front roller on the right-hand side.
11. Lift right hand end of the front roller enough to remove the drive belt.
12. Reverse the procedure to install the new drive belt.
13. After installation of the new drive belt, perform "Running Belt Tracking Adjustment" on page 71 and "Running Belt Tension Adjustment" on page 72.

Component/FRU Removal and Replacement

Bracket Motor Right Hand

Motor Tension

Drive Belt

(6) Bolt 1/4"-20 x 1 1/4" Long

(6) Bolt 1/4" Flat Washer

(6) Bolt 1/4" Lock Washer

(6) Bolt 1/4"- 20 Nut

Replacing/Adjusting the Elevation Actuator

Tools Required

- 2 - 9/16" Wrenches or Sockets
- 3/8" Wide Flat Blade Screwdriver
- Shipping blanket, carpet, or cardboard
- 1/4" Nut Driver or Flat Blade Screwdriver
- Grounding Wrist Strap
- Tape measure
- Service grease (Part Number 2097829-072)

Approximate time to complete: 1 hr 30 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Move the treadmill to an area large enough to lay it on its side.
5. Remove the (5) hood securing bolts.
6. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

7. Unplug the treadmill actuator harness.
The harness is a six-conductor connector attached to the actuator.
8. Cut the zip tie that holds the wire to the motor pan gusset.
9. With the shipping blanket in suitable spot, place a foot on the treadmill side channel, grasp the handrail; pull the treadmill towards you, and allow it to gently come to rest on its side.
10. With the 9/16" wrench and socket, remove the top and lower bolts from the actuator.
When removing the bolts from the actuator, watch out for nylon washers on each side of the actuator.
11. Prior to installing the new actuator, use the following procedure to preset the zero position.
 - a. Looking from the bottom (tube end), lightly turn the tube clockwise until it bottoms the threads or stops.


Component/FRU Removal and Replacement

- b. Ensure that the motor shaft does not turn.
 - c. Turn the tube counterclockwise 1.5 turns (540°) and continue turning it until the bolt holes in the tube align with holes as it is installed.
12. Reverse the procedure to install the actuator into the treadmill.
Be sure to use Service grease (quality, moly-based grease) on both the bolts and nylon washers retaining it.

NOTE:

Over tightening the lower bolt will lead to a popping or creaking noise when the treadmill is elevating.

13. After installing the elevation actuator, perform "Calibrating Elevation" on page 209.


Replacing the Speed Sensor from 2 Wire to 3 Wire

Tools required

- 1/4" Nut driver or flat blade Screwdriver
- Long Nose Pliers
- Grounding Wrist Strap

Approximate time to complete: 20 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the 2 wire connector from J10.


5. Feed wire connector thru hole on actuator support.

NOTE:

Remove two wire speed sensor and rubber grommet if it is a Quantum/TruTech configuration.

6. Slide 3 Wire Speed Sensor thru rubber grommet.


NOTE:

Depends on Motor type where Speed Sensor is routed and mounted.

7. Adjust the speed sensor to the magnet on the flywheel to approximately 1/4"

Kollmorgen

Quantum/TruTech

8. With the long nose pliers remove the DC Fan negative from J2.


Component/FRU Removal and Replacement

Caution: Use ESD precautions including wearing a grounding wrist strap which is connected to an exposed metal connection point on the treadmill chassis.

9. Reconnect the DC Fan negative to 3 Wire Speed Sensor Harness receptacle tab.

10. With long nose pliers, gently attach the 3 Wire Ground to J2 on the Smart Power Supply Board.


11. After replacing Speed Sensor, Perform Speed Calibration.

Replacing the Speed Sensor 3 Wire to 3 Wire

Tools required

- 1/4" Nut driver or flat blade screwdriver
 - Grounding Wrist Strap
- Approximate time to complete: 20 min
1. Turn the treadmill main power switch to the OFF position.
Unplug the treadmill power cord from the wall receptacle.

2. Allow the treadmill to discharge for 1 minute.
3. Remove existing 3 wire sensor.
4. Slide replacement 3 Wire Speed Sensor thru rubber grommet.


Caution: Use ESD precautions including wearing a grounding wrist strap which is connected to an exposed metal connection point on the treadmill chassis.

5. Adjust the speed sensor to the magnet on the flywheel to approximately 1/4"

Kollmorgen

Quantum/TruTech

6. After replacing Speed Sensor, Perform Speed Calibration.

Replacing the Smart Power Supply Relay Board

Tools required

- Wire Cutters or Utility Knife
- Grounding Wrist Strap
- 1/4" Nut Driver or Flat Blade screwdriver
- 5/16" Nut Driver
- 3/32" Flat Blade Electronics Screwdriver
- Grounding Wrist Strap

Approximate time to complete: 1 hr 15 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.

Component/FRU Removal and Replacement

4. Remove the (5) hood securing bolts.
5. Remove the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

6. Unplug the 4-conductor JP2 and 3-conductor JP1 plugs from the SPSR Board (smart power supply relay board).
7. Unplug the Cooling Fan J1-Red and J2-Black single conductor spades from the SPSR Board.
8. Unplug the E-Stop J4-Black and J5-Gray single-conductor spades from the SPSR Board.
9. Unplug the Stop Tether Switch J6-Black and J7-Gray single-conductor spades from the SPSR Board.
10. Unplug the Speed Sensor J10 two-conductor spade from the SPSR Board.
11. Unplug the Elevation Count Sensor J11-three-conductor spade from the SPSR Board.
12. Unplug the Drive Harness J15-ten-conductor spade from the SPSR Board.
13. Remove the RS232 harness J13 by unscrewing the (2) flat head screws retaining it to the SPSR Board.
14. Cut the Zip tie on USB harness J14 and unplug from the SPSR Board.
15. Remove the two control board bracket screws.
16. Remove the screw that holds the circuit board to the control board bracket.
17. Relocate the standoff and spacer transfer to the new Power supply board.
18. Remove the nylon nuts on the four corners of the SPSR Board.
19. Once the board has been removed, note the position of the voltage switch.
20. Set the switch on new board to the same voltage as the board just removed.

WARNING

Setting a different voltage on the new board will result in damage!!

21. Reverse the procedure to install the new SPSR Board.

22. After installation, see Chapter 9, "Calibration", to recalibrate the treadmill using the recalibration instructions for your treadmill model, and then check for proper operation.

Nylon Screw & Standoff

Control Board Bracket

Control Board Nut

Control Board Screw

Replacing the Circuit Breaker

Tools Required

- 1/4" Nut driver or flat blade screwdriver
- 3/32" Flat blade electronic screwdriver
- Grounding Wrist Strap

Approximate time to complete: 20 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (2) screws retaining the circuit breaker mounting plate.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

5. Gently pull the plate out, allowing the wires to feed through the hole.
6. Unplug the (2) wire connectors from the circuit breaker to be replaced.


When replacing more than one circuit breaker, perform the replacement for one breaker at a time to keep the wires on the correct value of circuit breaker.

7. Using a small electrical flat blade screwdriver, push in the holding clips on the sides of the breaker to allow it to slide out of the plate.

Component/FRU Removal and Replacement

The clips are part of the breaker housing and are located on the back side of the plate.

8. Install the new circuit breaker by sliding it into the hole until it snaps into place.
9. Reattach the wire connectors by sliding them on to the terminals.
10. Once all desired breakers have been replaced, reinstall the plate with the (2) screws.
11. After installation, check the unit for proper operation.

Relocating the Emergency Stop Button (ESB) Assy and Stop Tether Switch (STS) Assy

Tools Required

- 1/8" Allen Wrench
- 1/4" Flat Blade Screwdriver
- 2 Pair Long Nose Pliers
- Grounding Wrist Strap


Approximate time to complete: 1 hr 00 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. With the 1/8" Allen wrench, remove the (4) 10-32 screws (Fig. A) that secure the ESB and STS assembly to the fabricated handrail.

Fig. A

Fig. B

5. With the flat head screwdriver, lift and remove the ESB and STS assembly from the fabricated handrail (Fig. B).

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.


6. Completely remove the ESB and STS assembly and gently pull to expose the wires.
7. With two pairs of long nose pliers, gently separate the connectors on the ESB.
8. With a pair of pliers, remove the STS connector from the Stop Tether Switch Assy. (Fig. C).

Fig. C

Component/FRU Removal and Replacement

Fig. D

9. Reverse steps 1 through 8 to reinstall the ESB and STS assembly.
10. Remove the (5) hood-securing bolts.
11. Remove the hood and set it aside.
12. Swap J4 and J5 with J6 and J7 on the controller board.
With the pliers, cautiously remove the wires. Slowly wiggle the wires up and down while pulling. Avoid causing excessive flexing to the board (Fig. D).
13. Replace the hood and secure it with the original screws.

Replacing the Emergency Stop Button (ESB) Assy

Tools Required

- 1/8" Allen Wrench
 - 1/4" Flat Blade Screwdriver
 - 2 Pair Long Nose Pliers
 - Grounding Wrist Strap
- Approximate time to complete: 30 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. With a 1/8" Allen wrench, remove the (2) 10-32 screws (Fig. A) that secure the ESB assembly to the fabricated handrail.

Fig. A



Fig. B

5. With a flat head screwdriver, lift and remove the ESB assembly from the fabricated handrail (Fig. B).

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

6. Completely remove the ESB assembly and gently pull to expose the wires.
7. With two pairs of long nose pliers, gently separate the connectors on the ESB.
8. Reverse step 1 through step 7 to reinstall the ESB assembly.

NOTE:

When reconnecting the wires, they may be connected without regard to polarity.

Component/FRU Removal and Replacement

Replacing the Stop Tether Switch (STS) Assy

Tools Required

- 1/8" Allen Wrench
 - 1/4" Flat Blade screwdriver
 - Long Nose Pliers
 - Grounding Wrist Strap
- Approximate time to complete: 30 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow treadmill to discharge for 1 minute.
4. With a 1/8" Allen wrench, remove the (2) 10-32 screws (Fig. A) that secure the STS assembly to the fabricated handrail.

Fig. A



Fig. B

5. With a flat head screwdriver, lift and remove the STS assembly from the fabricated handrail (Fig. B).

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

6. Completely remove the STS assembly and gently pull to expose the wires.
7. With a pair of long nose pliers, remove the two connector wires from the outside terminals of the STS Assy. (Fig. C).

NOTE:

When reconnecting wires be sure to locate the connectors on the outer 2 terminals.


Fig. C

8. Reverse step 1 through step 7 to reinstall the STS assembly.

Replacing the Right or Left Handrail

Tools Required

- 1/8" Allen Wrench
 - 5/32" Allen Wrench
 - 5/16" Allen Wrench
 - Flat head screwdriver
 - Grounding Wrist Strap
- Approximate time to complete: 30 min

NOTE

When replacing all the handrails, remove the center handrail before proceeding with the left and right handrails. See "Replacing the Center Handrail" on page 189 for instructions.

1. Turn the treadmill's power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the ESB or STS assembly.

Component/FRU Removal and Replacement

NOTE:

(See "Replacing the Emergency Stop Button (ESB) Assy)" on page 184 or "Replacing the Stop Tether Switch (STS)" on page 186 for instructions.

5. Pay close attention not to damage the wires going through the handrail assembly.
6. With the 5/32" Allen Wrench, remove the four 1/4"-20 screws holding the center handrail to the handrail assembly.
7. Feed the wire harness through the handrail assembly.


NOTE:

Pay close attention not to damage the wires going through the handrail assembly.

- a. If replacing the right handrail, only disconnect the right side ESB (normal configuration) or STS (alternate configuration) wires to remove the handrail.
 - b. If replacing the left handrail, first disconnect the ESB (normal configuration) or STS (alternate configuration) on the right side and feed the wire harness through the center handrail. Then disconnect the STS (normal configuration) or ESB (alternate configuration) wires on the left side to remove the handrail.
8. With the flat head screwdriver, remove the two plastic caps covering the bolts on the handrail.
 9. With the 5/16" Allen Wrench, remove the two 3/8"-16 screws from the handrail.

This allows the handrail to pivot to a resting position beside the treadmill.
 10. Remove the remaining 3/8"-16 screws and nylon washer.

The handrail is ready to be replaced.
 11. Reverse steps 1 through step 10 to reinstall the replacement handrail.

Replacing the Center Handrail

Tools Required

- 1/8" Allen Wrench
- 5/32" Allen Wrench
- Grounding Wrist Strap

Approximate time to complete: 30 min

1. Turn the treadmill's power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the ESB or STS.

(See "Replacing the Emergency Stop Button (ESB) Assy)" on page 184 or "Replacing the Stop Tether Switch (STS)" on page 186 for instructions.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.


5. With the 5/32" Allen Wrench, remove the eight 1/4"-20 screws holding the center handrail to the handrail assembly.
6. Feed wire harness through center handrail.

Component/FRU Removal and Replacement

Pay close attention when removing as to not damage the wires running through the center handrail assembly.

The center handrail is ready to be replaced.

7. Reverse steps 1 through step 6 to install replacement center handrail.

Removing and Reinstalling the Handrails for Moving

Tools Required

- 1/4" Hex Socket or Flat Blade screwdriver
- 5/16" Allen Wrench
- Flat head screwdriver
- Grounding Wrist Strap

Approximate time to complete: 45 min

NOTE:


This task is required only when you need to pass through a door opening less than the standard 36". It will be necessary to remove the handrail assembly if the door is not capable of being opened fully parallel to the door opening and clearance is less than 35 1/2".

1. Turn the treadmill's power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (5) hood securing bolts.
5. Removed the hood and set it aside.

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis.

6. At the Smart Power Supply Board, disconnect the two wire harnesses from the J4 and J5 connectors (ESB) and the J6 and J7 connectors (STS).
7. Slide the ESB and STS wires through the side channel of the motor pan.
Pay close attention not to damage the wires going through the handrail assembly.
8. With the flat head screwdriver, remove the two plastic caps on each side which are covering the bolts on the handrails.
9. With the 5/16" Allen Wrench, remove the two 3/8"-16 screws from each handrail.
This allows the handrail to pivot to a resting position beside the treadmill.
10. Remove the remaining 3/8"-16 screws and nylon washers.
The handrail is ready to be lifted off as an assembly.
11. Reverse step 1 through step 10 to reinstall the replacement handrail.

Replacing the Main Power Switch

Tools Required

- 1/4" Nut driver or flat blade screwdriver
- 3/32" Flat blade electronic screwdriver
- Grounding Wrist Strap

Approximate time to complete: 25 min

1. Turn the treadmill main power switch to the OFF position.
2. Unplug treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.

Component/FRU Removal and Replacement

CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis


4. Remove the two screws retaining the main power switch mounting plate.
5. Gently pull the plate out, allowing the wires to feed through the hole.

NOTE:

If Switch plate is damaged during switch replacement see 5a or 5b below for proper GE FRU part numbers to use for identified Serial Number(s).


2097829-039 for GEDC-6608 and earlier production and 2097829-124 for GEDC-6609 and later production.

- a. Installation base prior to GEDC-6608 use GE FRU Part Number 2097829-039 (MOTOR PAN HARDWARE KIT)


- b. Forward production after GEDC-6608 use GE FRU Part Number 2097829-124 (MTR PAN HARDWARE KIT MC PRODUCT & ALT SCHURTER PLATE).

5b-1 – Schurter Cover Plate + GE Part Number 2097829-036 (SWITCH TA35-C324-L200C0) if switch is needed to be replaced.

5b-2 – Mechanical Products Cover Plate + GE Part Number 2097829-123 (SWITCH MECHANICAL PRODUCT 2421-305-119-2000) if switch is needed to be replaced.

6. Unplug the four (4) wire connectors from the main power switch to be replaced.

CAUTION:


When reconnecting the main power switch, insert the brown wire on top and black wire on bottom, as shown in the picture in step 8.

7. With a small electrical flat blade screwdriver, push the holding clips on the sides of the main power switch to allow it to slide out of the plate.

The clips are part of the power switch housing and located on the back side of the plate.

Component/FRU Removal and Replacement

8. Reverse step 1 through step 7 to reinstall main power switch while being sure to follow the caution in step 6.

Adjusting the Run Belt Tracking and Tension

Tools Required

- 1/4" Allen wrench

Approximate time to complete: 25 min

Adjusting belt tracking

Your running belt should remain centered, although a slight amount of movement to the left or right during use is normal. Do not allow the running belt to travel all the way to either side.

To adjust the belt tracking, do the following:

1. Turn the treadmill's power switch to the ON position.
2. Increase the speed to 3 mph (4.8 km/h)
3. Observe the left side of the running belt as it travels over the rear roller.

If the belt runs to the right side of the roller, turn the right bolt 1/8 turn clockwise and turn the left bolt 1/16 turn counterclockwise.

NOTE:


When tightening one side of the belt, always loosen the opposite side one-half as much. This procedure provides finer control, with a smaller impact on the belt tension.

4. Check the belt after 2 minutes, with the treadmill running at approximately 7 mph (11.3 km/h).

If the belt does not correct itself, continue with slight turns until the belt is in the center of the rear roller. If the belt runs toward the left side of the roller, reverse the adjustments.

NOTE:

Uneven floors accelerate belt misalignment. This situation may require more frequent adjustment to prevent belt damage.

Adjusting belt tension

The running belt may stretch and loosen with regular use. This looseness is noticeable when the belt tends to hesitate or stick.

To adjust the belt's tension, do the following:

1. Turn the treadmill's power switch to the ON position.
2. Start the treadmill and increase to 1 mph (1.6 km/h).
3. Start walking on the treadmill, grabbing the side handrail and applying pressure with your foot to create resistance on the running belt.
4. If the running belt hesitates or slips on the front drive roller, tighten both tension bolts 1/2 turn (clockwise).
5. Repeat steps 2 through step 4 until the running belt stops slipping.

NOTE:

When the running belt is too tight, the edge of belt will curl, causing premature running belt failure.

Replacing the Power cord


Tools Required

- 7/16" Open End Wrench or Socket

Approximate time to complete: 10 min

1. Turn the treadmill's power switch to the OFF position.
2. Unplug the treadmill from the wall receptacle.
3. Allow treadmill to discharge for 1 minute.
4. With the 7/16" wrench, remove the 1/4"-20 bolt, allowing the clamp to be removed from the front of the treadmill.
5. Unplug the C-19 appliance power cord from the front of the treadmill.
6. Reverse step 1 through step 5 to reinstall new power cord.

Component/FRU Removal and Replacement

Replacing the CE filter

Tools Required

- 1/4" Nut driver or flat blade screwdriver
- Grounding Wrist Strap
- Pliers or adjustable wrench

Approximate time to complete: 30 min


1. Turn treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (5) hood-securing bolts.
5. Remove the hood and set it aside.

CAUTION:


Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis. Put on the grounding wrist strap and ground it to the power cord ground screw.

6. With pliers or an adjustable wrench, hold the 11/32" nuts while removing the two (2) #8-32 screws with a flathead screwdriver.
7. Remove the quick disconnect wire connection from the inline filter.

Note the wire's location for replacement.

8. Reverse step 1 through step 7 to reinstall replacement CE Filter.

Replacing the Relay

Tools Required

- 1/4" Nut driver or flat blade screwdriver
- Grounding Wrist Strap
- Pliers or adjustable wrench

Approximate time to complete: 20 min

1. Turn treadmill main power switch to the OFF position.
2. Unplug the treadmill plug from the wall receptacle.
3. Allow the treadmill to discharge for 1 minute.
4. Remove the (5) hood-securing bolts.
5. Remove the hood and set it aside.


CAUTION:

Use ESD precautions, including wearing a grounding wrist strap that is connected to an exposed metal connection point on the treadmill chassis. Put on the grounding wrist strap and ground it to the power cord ground screw.

6. With pliers remove the quick disconnect wires connection from the relay.

Note the wire's location for replacement.

Component/FRU Removal and Replacement

7. Reverse step 1 through step 6 to reinstall replacement relay.

9

Calibration

Calibration is necessary to ensure your T2100-ST series treadmill is functioning properly. This chapter describes how to:

- install the T2100-ST series calibration software,
- install the serial drivers,
- connect the calibration software to the T2100-ST series treadmill,
- calibrate the treadmill's speed,
- calibrate the treadmill's elevation, and
- understand the calibration software's tabs.

Installing the T2100-ST Series Calibration Software

This section explains how to install the T2100-ST Series Calibration Software on the field service engineer's PC in preparation for calibrating the T2100-ST Series treadmill. Compatible PC operating systems are Windows 7, 8, 8.1, and 10 (both 32 bit and 64 bit).

NOTE:

This section can be skipped if the T2100-ST Series Calibration Software is already installed on the PC being used. If already installed, proceed to "Installing the Serial Drivers" on page 201.

NOTE:

Do not plug any hardware including treadmill cable into the PC unless directed to do so.

1. Download the T2100-ST Series Calibration software using the following procedure:
 - a. Go to the GE Service Support Central site;
http://supportcentral.ge.com/products/sup_products.asp?prod_id=20335

Calibration

NOTE

A Single Sign On (SSO) user ID and password are required to log on to the site.

- b. In the **Product Support** box, select **Carts-Stress Software** under the **Tools** heading.
GE Libraries opens to the **Carts-Stress** folder.
- c. Select **Software Library > T2100 ST1 and ST2 Treadmill** and download the correct version of the calibration software.

There are two versions of calibration software:

- T2100-ST_Series_Calibration_Software_1_0_1_5.zip This version of the software is for calibrating firmware version 1.57. If the T2100-ST treadmill Firmware has been upgraded to V1.66 this software is no longer compatible.
- T2100 ST1 and ST2 Calibration Software _1_0_1_28.zip This version of the software is for calibrating firmware version 1.59+ thru 1.18.12+ & 1.19.109+

For more information, refer to *T2100 ST1 and ST2 Calibration Software Information.docx*, also found at this location.

NOTE


Service personnel can also call Technical Support to receive the calibration software via email.

2. To install the T2100-ST Series calibration software, extract the zip file and double-click **setup.exe** as located in the previously chosen directory.

NOTE:


The following instructions are for a PC running Windows 7. The basic instruction flow should be the same for Windows 8 and 10, although the screens may vary. If necessary, refer to Appendix 0, "Reference Documents", for a list of reference documents that contain helpful OS-specific guidance.

The **Application Install – Security Window** opens.

3. Click **Install**.

With a successful software installation, the software starts automatically.

4. In preparation for the next section, close the **T2100-ST Series Calibration Software** window.

Installing the Serial Drivers

This section describes how to install the USB Serial Drivers on the field service engineer's PC in preparation for running the Calibration Software. Compatible PC OSes are Windows 7, 8, 8.1, and 10 (both 32 bit and 64 bit).

NOTE:

This section can be skipped if the USB Serial Drivers are already installed on the PC being used. If already installed, proceed to "Software Interface Instruction T2100-ST Series Calibration" on page 205.

Calibration

NOTE:

Do not plug any hardware including the treadmill USB cable into the PC until this step is completed.

1. Locate the flash drive (p/n 317-160-313) shipped with the T2100-ST series treadmill and insert into a USB port of the PC.
2. Browse the flash drive and locate the file named **CDM21218_Setup.exe** and double click on it to start the install.

NOTE:


Do not download the driver from the FTDI website.

NOTE:

The following instructions are for a PC running Windows 7. The basic instruction flow should be the same for Windows 8 and 10, although the screens may vary. If necessary, refer to Appendix O, "Reference Documents", for a list of reference documents that contain helpful OS-specific guidance.

3. If prompted 'Do you want to allow the following program to make changes to this computer', verify the program name, **CDM21218_Setup.exe**, and publisher, **Future Technology Devices International, Ltd**, and click **Yes**.

The **Driver Extract** window opens.


4. Click **Extract**.

The **Device Driver Installation Wizard** opens.

5. Click **Next**.
The **License Agreement** window opens.

6. Select **I accept this agreement** and click **Next**.
An installation progress window opens.

Calibration

The installation process should take less than a minute.

The **Device Driver Installation Wizard** window opens, confirming a successful installation.

7. To close the window, click **Finish**.

Software Interface Instruction T2100-ST Series Calibration

Tools Required

- Calibration Program (T2100-ST Series Calibration Software)
- Chalk or some other temporary marking device
- Computer or laptop (PC)
- Serial or USB Communication Cable

NOTE:

This task needs to be performed prior to speed or elevation calibration.


1. Connect the treadmill to AC power.
2. Connect the treadmill and PC.
 - If using USB cable, connect a USB cable (GE part number 2044095-001) between the PC and the T2100-ST series treadmill and turn on the treadmill.
 - If using RS232 serial cable:
 - a. Connect one end of the RS232 serial cable (GE part number 2097829-109) to the RS232 COM port located on the lower front panel of the T2100-ST series treadmill.
 - b. Connect the other end of the RS232 serial cable (GE part number 2097829-109) to a COM port on the PC.
 - c. Turn on the Treadmill.
3. Open T2100-ST Series Calibration Software by browsing for the file named **T2100-ST_Series_Calibration_Software.application** and then double clicking on it.

NOTE:

This file should be found in the directory in which the file **ST_Series_Calibration_Software.Zip** was unzipped.

4. In the lower left of the window, select the appropriate COM and press **Open**.

For example:

5. When connected, the numerical count to the right of **Sent:** will increase in value and detect the drive/motor configuration:

ReScan 1.0.1.27 Sent:2057 M/C Type : Quantum

M/C Type can also be Kollmorgen for older configurations.

This is the indication that the communication has been established with the treadmill.

Calibration

Calibrating Speed


Tools Required

- Calibration Program (T2100-ST Series Calibration Software)
- Stopwatch (handheld or software app) or watch with a second hand
- Chalk or some other temporary marking device
- Computer or laptop
- Serial or USB Communication Cable

Approximate time to complete 20 min

1. If not already using it, open the T2100-ST Series Calibration Software.
2. In the lower left of the window, select the appropriate COM and press **Open**.

For example:

When connected, the numerical count to the right of **Sent:** will increase in value:

ReScan 1.0.1.27 Sent:2057

This is the indication that the communication has been established with the treadmill.

3. Select the **Advanced Control** tab.
Select the **Set PC Time**.

NOTE

When replacing the 4th edition Smart Power Supply Board Set PC Time is required. When updating a 3rd edition Smart Power Supply Board this feature is not available.

4. Select the **Speed Cal** tab.
5. Remove all tools from the running belt and stand clear of the running surface.
6. Select **Reset Treadmill**.

The treadmill will reset and return to the parked position.

NOTE

After the down relay on the Smart Power Supply board cycles (**TWO AUDIBLE CLICKS**), a 20 second timer starts during which the following step (start speed calibration) must be initiated. If this time window is missed, this step must be repeated.

7. Select **Start Speed Cal.**

You will be promoted with "Treadmill Entering Calibration".

8. Click **OK**.

The running belt will pause for approximately 8 seconds and slowly increase speed, establishing multiple calibration points. This procedure may take approximately 10 to 15 seconds per calibration point, with the overall time of the speed calibration procedure taking up to 5 minutes to complete.

When calibration routing is complete the treadmill will stop and, if successful, the Speed Calibration table will have automatically updated all of the values in the **Measured Speed** column (for example see table below). At this point running belt calibration is complete.

PWM Output	Measured Speed
PWM Output 1%	0.151
PWM Output 5%	0.774
PWM Output 10%	1.552
PWM Output 20%	3.112
PWM Output 30%	4.679
PWM Output 40%	6.252
PWM Output 50%	7.834
PWM Output 60%	9.426
PWM Output 70%	11.026
PWM Output 80%	12.641
PWM Output 85%	13.252
PWM Output 87.5%	13.266
PWM Output 90%	13.285
PWM Output 93.5%	13.298

NOTE:

If the MEASURED SPEED field are displaying "0" post calibration, this indicated that the calibration failed, refer to the flow chart titled **Smart Power Supply Error Code 1 Flow Chart 1J "Bad Calibration Error (1)"** on page 100.

After replacement of the DRIVE BOARD or MOTOR a new SPEED calibration is required.

9. Select the **Basic Control** tab, select **Start Belt**, and then press **▲** one time.

The running belt will start at 0.1 mph.

NOTE

You have two ways to change speed:

- selecting **▲** or **▼** and adjusting to the desired speed, or

Calibration

- double-clicking on the **Set Speed** and typing in the desired speed

Set Speed

0.10

10. To validate the belt speed calibration, do the following:
 - a. Verify a point on the run belt that can be seen plainly with each full revolution; this could be the seam or a chalk line you make on the belt.
 - b. Start the treadmill and adjust the speed to 1 mph/1.6 km/h.
 - c. Start the stopwatch when the mark on the belt passes an established point on the treadmills side rail.
 - d. Count the number of passes the mark makes in a 60 second period. The count should be 8 passes; this is within the 2% speed accuracy.
 - e. Set the speed to 5 mph/8.0 km/h.
 - f. Count the passes the belt makes but this time in a 15 second period. The count should be 10 passes; this is within the 2% speed accuracy.
 - g. Do one of the following, depending on the speed accuracy:
 - If treadmill achieved 2% accuracy, select **Stop Belt**.
Speed calibration is complete.
Proceed to "Calibrating Elevation" on page 209 if a Smart Power supply or elevation actuator was replaced.
 - If the treadmill does not achieve the 2% speed accuracy, this can be adjusted. When needing to further adjust the speed accuracy, this can be achieved by changing the Motor Controller Speed Offset and External Sensor Speed Offset.

Proceed to step 11.

11. Select **Stop Belt** and return to the **Speed Cal Tab**.
12. If the actual running belt is running slower than target speed, double-click the **Motor Controller Speed Offset** and decrease the percentage between 95 to 99.99%.
13. If the actual running belt is running faster than target speed, double-click the **Motor Controller Speed Offset** and increase the percentage between 100.01 to 104.99%.
14. Adjust the **External Sensor Speed Offset** to match the changed percentage of the **Motor Controller Speed Offset**.

This is critical to ensure the treadmill does not shut down due to incorrectly reported speed verification (refer to Smart

Power Supply Error Code 3 through 7 in Chapter **Error!**
Reference source not found., **"Error! Reference source not found."**.

15. Repeat steps 6 through 10 to confirm proper adjustment of the actual running belt speed.

Calibrating Elevation

NOTE:

Calibrating Elevation is required if a Smart Power Supply or Elevation Actuator was replaced, or if reported elevation is not achieving desired accuracy.


Tools Required

- Calibration Program (T2100-ST Series Calibration Software)
- Tape Measure
- 1/4" Nut Driver or Flat Headed screwdriver
- T15 Torx screwdriver
- Phillips Head screwdriver
- Needle Nose Pliers
- Computer or laptop
- Serial or USB Communication Cable

Approximate time to complete 30 min

1. If not already using it, open the T2100-ST Series Calibration Software.
2. In the lower left of the window, select the appropriate COM and press **Open**.

For example:

When connected, the numerical count to the right of **Sent:** will increase in value:

ReScan 1.0.1.4 Sent:6091

This is the indication that the communication has been established with the treadmill.

3. Select the **Grade Cal** tab.
4. Select **Test 0.5%** to elevate the treadmill 0.5% and then cycle the treadmill power.
By resetting the power, the treadmill will automatically find the lower limit position.
5. Using a tape measure, verify the 0% elevation.

Calibration

Refer to the "Elevation Chart" on page 212.

- a. If the zero position requires adjustment, remove actuator top cover.

VonWeise Actuator

MSI Actuator

**NOTE:**

The system can have either VonWeise or MSI Actuator they only difference is the zero position is adjusted.

- b. Do one of the following:
 - If the position is low, adjust the upper cam (control for 0% elevation limit) with needle nose pliers clockwise into the limit switch "see note".
 - If position is high, adjust the upper cam (control for 0% elevation limit) with needle nose pliers counterclockwise away from limit switch "see note".

VonWeise Actuator

MSI Actuator

**NOTE:**

When adjusting upper cam on VonWeise Actuator, an audible click will be heard. Make adjustments in small increments, preferably a single click at time.

On MSI Actuator loosen screw to relieve tension then position upper cam. This is a friction fit be careful to move in small increments, tighten screw when complete.

- c. Repeat step 4 to 5 to reestablish 0% elevation.
6. To ensure the upper limit switch is not engaged, change the 25% default value of **481** on the **Grade Cal** tab to a new value of **500**.
7. Select **Test 25%** and verify with a tape measure that the value is 18.0" to 18.1".

The lower cam should almost be in contact with the limits switch.


If the **Upper Elevation Limit** is reached prior to reaching the target, the treadmill will automatically reset to the 0% elevation.

- a. If this happens on VonWeise Actuator, use a large (.375" wide) flat blade screwdriver to rotate the lower switch cam counterclockwise the proper amount to assure the treadmill will not elevate past 25%. On MSI Actuator loosen screw to relieve tension then position lower cam. This is a friction fit be careful to move in small increments, tighten screw when complete.

VonWeise Actuator

MSI Actuator

- b. If the cam was turned to adjust the lower switch (control for 25% elevation limit); the upper cam (control for 0% elevation limit) will need to be reset to zero to allow the treadmill to park in the correct position. Repeat step 4 to 5 to reestablish 0% elevation.

NOTE:

When an adjustment on VonWeise Actuator is made to the lower cam, the upper cam moves in tandem.

On MSI Actuator when an adjustment is made to the lower cam, the upper cam does not move in tandem.

THE LOWER SWITCH adjusts the UPPER LIMIT and the UPPER SWITCH adjusts the LOWER LIMIT. If the UPPER LIMIT has been adjusted, then the LOWER LIMIT will also need to be RE-ADJUSTED and vice versa

8. When the desired results are achieved change the **Test 25%** value back to **481**.
9. Select **Test 0.5%**.
10. Using the tape measure, verify the elevation height Refer to the "Elevation Chart" on page 212.
11. If incorrect, compensate by making small incremental adjustments of 1 to 5 counts in **Elevation Count** column next to the elevation grade being tested by typing the new value and selecting **OK** or pressing **Enter**.

Grade	Inches	MM	Elevation Count	Test
0%	2.64"	67mm	0	Home
0.5%	3.35"	85mm	15	Test 0.5%
5%	5.92"	150.4mm	93	Test 5%
10%	9.18"	233.2mm	134	Test 10%
15%	12.45"	316.2mm	296	Test 15%
20%	15.12"	384mm	389	Test 20%
25%	17.88"	454.2mm	481	Test 25%

12. Select **Test 0.5%** and recheck the measurement.


Calibration

NOTE:

When selecting the grade being tested, the treadmill will change elevation.


13. Repeat steps 9 through 12 for 5%, 10%, 15%, 20%, and 25%.
14. Cycle the power to return the treadmill to 0% elevation.
15. If the actuator cap was removed, replace the cap.

Elevation Chart


Elevation Actuator

VonWeise Actuator

MSI Actuator

Calibration Software Tabs

Basic Control Tab

Calibration

Start Running Belt: Select **Start Belt**

Speed Control: Press ▲ one time. The running belt will start at 0.1 mph. You have two ways to change speed:

- Select ▲ or ▼ and adjust to the desired speed
- Double-click on **Set Speed** and type in the desired speed

Set Speed

0.10

Stop Running Belt: Select **Stop Belt**. The running belt will go to and maintain 0.0 mph for about 3 minutes before the belt starts to freewheel.

MPH Versus KPH: The software is configured to run in MPH as a standard configuration. When it is desired to run in KPH, check the **KPH** box. The next time you select **Start Belt**, the treadmill will function in KPH.

Grade Control: Press ▲ one time, the elevation will increase to .5%. You have two ways to change grade:

- Select ▲ or ▼ to adjust elevation by .5% per click.
- Double click on **Set Grade** and type in the desired grade %.

Set Grade

0.50

COM port: Interface COM port with host system.

Advanced Control Tab

T2100 - ST Series 4th Edition Calibration Software - Target Version 1.18...
Basic Control **Advanced Control** Speed Cal Grade Cal Error Code System Log Quantum

Reset Treadmill

Load Defaults

Load Pre-Calibration

Runtime : 3:01:00

Firmware Version : 1.19.103

Real Time Clock : 7/27/2023 10:52:07 AM

Set PC Time

Disconnect COM1 ReScan 1.0.1.27 Sent:31845 M/C Type: Quantum

Reset Treadmill: Select **Reset Treadmill**. The treadmill will perform a software reset: it will return to 0.0% elevation and reset the drive.

Load Defaults: Select **Load Defaults** to reset the speed calibration to 0.000.

NOTE:

The treadmill needs speed calibration to start the running belt.

Load Pre-Calibration: Select **Load Pre-Calibration** to upload calibration values.

NOTE:

Only to be performed when calibrating speeds fails, to assist in troubleshooting.

Set PC Time: The FGLF0495-3 Smart Power Supply is equipped with a real time clock. The time and date are required to be updated by selecting Set PC Time to match the computer system. This will assist when troubleshooting is required.

NOTE:

The FGLF0495-1 Smart Power Supply is not equipped with real time clock.

Runtime: Information indicating treadmill total runtime.

Firmware Version: Software version loaded in the smart power supply board.

Treadmill running firmware Version 1.XX does not have the ability to Set PC Time displays as follows:

T2100 - ST Series 4th Edition Calibration Software - Target Version 1.12 +
Basic Control Advanced Control Speed Cal Grade Cal Error Code System Log

Reset Treadmill**Load Defaults****Load Pre-Calibration**

Runtime : 2 11:30
Firmware Version : 1 1900
Real Time Clock : 5/17/2021 1 35 12 PM

Set PC Time

Disconnect COM13 + ReScan 1 0 1 23 Sent 934 M-C Type Quantum

Calibration

Speed Cal Tab

T2100 - ST Series 4th Edition Calibration Software - Target Version 1.18...

Basic Control Advanced Control **Speed Cal** Grade Cal Error Code System Log Quantum

Reset Treadmill

Start Speed Cal

No Calibration Error

Motor Controller Output : 0.0%

Measured Speed Motor Controller : 0.000 mph

Motor Controller Speed Offset 101.00%

Measured Speed External : 0.000 mph

External Sensor Speed Offset 101.00%

PWM Output	Measured Speed
PWM Output 1%	0.100
PWM Output 5%	0.931
PWM Output 10%	1.853
PWM Output 20%	3.699
PWM Output 30%	5.545
PWM Output 40%	7.389
PWM Output 50%	9.236
PWM Output 60%	11.082
PWM Output 70%	12.886
PWM Output 80%	14.773
PWM Output 85%	15.597
PWM Output 87.5%	15.947
PWM Output 90%	16.197
PWM Output 93.5%	16.447

Disconnect COM1 • ReScan 1.0.1.27 Smb:31851 M/C Type : Quantum

Reset Treadmill: Select **Reset Treadmill**. The treadmill will perform a software reset: it will return to 0.0% elevation and reset the drive.

Start Speed Calibration: Select **Start Speed Cal** to initialize the speed calibration routine.

No Calibration Error: Information regarding calibration error.

Motor Controller Output: When the treadmill speed is running, this displays the motor controller % of speed output.

Measured Speed Motor Controller: Information from the motor controller speed sensor.

Motor Controller Speed Offset: The ability to fine tune speed output when outside the 2% speed accuracy.

Measure Speed External: External speed sensor is in place to detect speed abnormalities.

External Sensor Speed Offset: The ability to fine tune the speed verification to the Motor controller speed sensor.


PWM Output: Calibration information reported in % to measured speed.

Grade Cal Tab

T2100 - ST Series 4th Edition Calibration Software - Target Version 1.18...

Basic Control Advanced Control Speed Cal **Grade Cal** Error Code System Log Quantum

Grade	Inches	MM	Elevation Count	Test
0%	2.64"	67mm	0	Home
0.5%	3.35"	85mm	19	Test 0.5%
5%	5.92"	150.4mm	93	Test 5%
10%	9.18"	233.2mm	194	Test 10%
15%	12.45"	316.2mm	298	Test 15%
20%	15.12"	384mm	389	Test 20%
25%	17.88"	454.2mm	481	Test 25%

Disconnect COM1 • ReScan 1.0.1.27 Sent:31858 M/C Type: Quantum

Elevation Count: The ability to fine tune the actual elevation target height.

Test: The ability to test each elevation calibration point.

Error Code Tab

T2100 - ST Series 4th Edition Calibration Software - Target Version 1.18...

Basic Control Advanced Control Speed Cal Grade Cal **Error Code** System Log Quantum

Error Code : None

Error Type	Error Count
Bad Calibration Error (1)	0
Elevation Error (2)	0
Zero Speed Motor Controller Error (3)	0
Over Speed Motor Controller Error (4)	0
Zero Speed External Sensor Error (5)	0
Over Speed External Sensor Error (6)	0
Speed Compare Error (7)	0
Motor Controller Fault Signal (8)	0
Belt Start Reject Error (9)	0
CDSE Count	0
Reset Count	0
Pull Tether Count	0
CPU Oscillator Trap Count	0
CPU Address Trap Count	0

Clear Error Counts

Stop Reason : 0021
Found Zero Elevation after Powerup


Disconnect COM1 • ReScan 1.0.1.27 Sent:31866 M/C Type: Quantum

Error Code: The ability to view error codes that have accrued during operation for the treadmill.

Clear Error Counts: The ability reset the error code counter to zero.

Calibration

System Log Tab

Set Log File Location: The ability to record treadmill communication for trouble shooting.

Log Events to File: When checked this will create a TXT file for trouble shooting. When not selected, it only displays information in the window.

Quantum Tab

MC Parameter	MC Value
MC Version	1280
MC UVW Fail Count	0
MC PWM	0
MC RPM	0
MC Current	0
MC Bus Volts	320
MC Flag	0
MC State	36
MC Fault	0
Max Est Accel	0
Max MC Accel	0

Disconnect COM1 ReScan 1.0.1.27 Sent:31883 M/C Type: Quantum

MC Parameter: The motor control variable measured by the calibration software.

MC Value: Value of the measured variable in *MC Parameter*.

10

Functional / Post Repair / Preventative Maintenance Checkout Procedures

This chapter identifies the tools and checkout procedures that need to be performed to ensure the T2100-ST Series treadmills and individual components are functioning properly and safe to use. These procedures are performed after completing any repair, exchange, or preventative maintenance procedure.

Tools Required

All or some of the following tools may be required to perform the various checkout procedures. For a list of which tools are required for each procedure, refer to "Procedural Lists" on page 225.

1. Applicable service and/or operator manual, as needed.
2. Standard hand tools.
3. Anti-static wrist strap.
4. Leakage tester.
5. Multimeter.
6. T2100-ST Series Calibration software.

Checkout Procedures

The following topics detail the individual checkout procedures. For a list of which procedures to perform after completing each repair or preventative maintenance procedure, refer to "Procedural Lists" on page 225.

Functional / Post Repair / Preventative Maintenance Checkout Procedures

Visual/Functional Checks

1. Perform the following inspections for visual signs of excess wear or damage.
 - a. Check the device case for cracks or any other visual damage.
 - b. Check plugs, cords, cable insulation jacket, and connectors for bent prongs or pins.
 - c. Check handles for wear and or damage.
 - d. Check walking belt for wear and or damage.

Expected Results: No visual damage and wear detected

Pass/Fail

Operational Checks

2. Inspect the following to ensure they are operational.

- a. Power Switch
 - o Connect T2100-ST Series to external power.
 - o Power on the device

Expected Results: System powers on

Pass/Fail

- b. Emergency Stop Switch
 - o Power on T2100-ST Series
 - o Adjust speed manually by controlling the unit.
 - o Press the Emergency Stop Button (ESB).
 - o To reset, turn the button 1/4-turn in clockwise direction or pull up.

Expected Results: Unit will coast to a stop; running belt will move freely.

Pass/Fail

- c. Stop Tether Switch (STS)
 - o Power on T2100-ST Series
 - o Adjust speed manually by controlling the unit.
 - o Pull the STS Clip
 - o To reset, replace the Clip or Magnet Tether in System.

Expected Results: Unit will have a controlled deceleration to a stop; running belt will not move freely.

Pass/Fail

Speed Calibration

3. Perform the following speed calibration and verification:
 - a. Power on T2100-ST Series treadmill.
 - b. Perform the "Calibrating Speed" procedure on page 206.

Expected Results: The system speed is successfully calibrated without errors.

Pass/Fail

Elevation Calibration

4. Perform the following elevation calibration and verification:
 - a. Power on T2100-ST Series.
 - b. Perform the "Calibrating Elevation" procedure on page 209.

Expected Results: The system elevation is successfully calibrated without errors.

Pass/Fail

Test Plug Check

5. Perform the following tests.
 - a. Turn off and remove the controlling device from back of treadmill.
 - b. Plug the Test Plug (secured to the frame) into the RS232 port.
 - c. Perform the "Test Plug Procedure" on page 50.
 - d. When the unit passes the test, turn it off and reconnect to controlling device.

Expected Results: The treadmill will elevate to 20% with a cycled speed of 10 mph, and then return to 0% and 0.0 mph.

Pass/Fail

Walking Belt Tension / Tracking

6. Perform the following to check the walking belt tension and tracking.
 - a. Power on the T2100-ST Series treadmill.
 - b. Adjust the speed manually by controlling the unit.
 - c. Perform the procedure "Running Belt Tension Adjustment" on page 53.
 - d. Perform the procedure "Secure Cables
 - e. Tie down cables to ensure they do not get caught in the wheels or the elevation racks.
 - f. Running Belt Tracking Adjustment" on page 53.

Expected Results: Belt should not slip and the rear roller should have no alignment tracking problem.

Pass/Fail

Walking Belt Speed Verification

7. Perform the following to check the walking belt speed.
 - a. Power on T2100-ST Series
 - b. Adjust speed manually by controlling the unit.
 - c. Set speed to 5.0 mph/8.0 km/h.

Expected Results: 10 revolutions of the running belt in 15 seconds.

Pass/Fail

Functional / Post Repair / Preventative Maintenance Checkout Procedures

Host Communication

8. Perform the following to test the host device's communication with the treadmill.
 - a. Attach the controller cable for the host system to the T2100-ST Series treadmill.
 - b. Power on the T2100-ST Series treadmill.
 - c. Power on the host system.
 - d. Adjust the treadmills speed and elevation from the host system.
Expected Results: Elevation and speed can be adjusted with the controlling device. **Pass/Fail**
 - e. Stop the treadmill from the host system.
Expected Results: Unit can be stopped with the controlling device **Pass/Fail**

Electrical Safety Checks

9. Perform the following to test the electrical safety of the treadmill.
 - a. Conduct the tests described in "Electrical Safety Tests" on page 49 and "Conducting Leakage Tests" on page 51.
 - b. Report the pass/fail status for each step performed.
 - c. Report the measured values.**Expected Results:** All electrical safety tests pass. **Pass/Fail**

Earth Leakage Current to Ground

					limits
1.	Forward Polarity (L1, L2)	NC	_____ μ A	Pass/Fail	500 μ A
2.	Neutral (L2) open, Forward Neutral Polarity	SFC	_____ μ A	Pass/Fail	500 μ A
3.	Ground open, Forward Polarity	SFC	_____ μ A	Pass/Fail	500 μ A
4.	Ground open, Reverse Polarity	SFC	_____ μ A	Pass/Fail	500 μ A
5.	Neutral (L2) open, Reverse Polarity	SFC	_____ μ A	Pass/Fail	500 μ A
6.	Reverse Polarity (L2, L1)	NC	_____ μ A	Pass/Fail	500 μ A

Ground Continuity

					Resistance
1.	AC mains power cord ground prong to exposed metal surface (ground lug)	N/A	_____ m Ω	Pass/Fail	100m Ω

Preventative Maintenance Visual/Functional Checks

10. Perform the following inspections for visual signs of excess wear or damage.
 - a. Check the device case for cracks or any other visual damage.
 - b. Check plugs, cords, cable insulation jacket, and connectors for bent prongs or pins.
 - c. Check handles for wear and or damage.
 - d. Press the Emergency Stop Button (ESB) and the treadmill will coast to a stop. To reset the ESB, turn the button 1/4-turn in clockwise direction or pull up. The treadmill will return to 0.0% elevation.
 - e. Pull the STS Clip. The treadmill will have a controlled stop. To reset, replace the Clip or Magnet Tether to the original position on the switch.
 - f. The treadmill will return to 0.0% elevation.
 - g. Evaluate the condition of the deck and running belt.
 - h. Check running belt tension and correct as required.
 - i. Observe running belt tracking and correct as required.
 - j. Check the drive belt tension and correct as required.
 - k. Clean and lubricate the treadmill elevation screw.

If the treadmill hood has been opened, check that the interior is free of excessive dust build-up, if necessary; use commercially available compressed air to clean following manufacturer's instructions.

Expected Results: No visual damage, no wear detected, maintenance checks / procedures performed successfully.

Pass/Fail

Procedural Lists

The following tables identify the tools, checkout procedures, and safety electrical checks required after completing each repair or preventative maintenance procedure. To use these tables, find the applicable procedure in the first column, then locate the **Tools**, **Functional Checkout Procedures**, and **Electrical Safety Checks** associated with that procedure. Use the numbers in those columns to locate the corresponding details earlier in this chapter.

For a list of the tools, refer to "Tools Required" on page 221.

For the checkout procedure details, refer to "Checkout Procedures" on page 221.

For the electrical safety check details, refer to "Electrical Safety Checks" on page 224.

Functional / Post Repair / Preventative Maintenance Checkout Procedures

FRU Repairs / Exchange

Description	Tools	Functional Checkout Procedure	Electrical Safety Checks
DC Motor Assembly Kollmorgan 2097829-020 and associated hardware.	1, 2, 3, 4, 6	1, 3, 8	9
PARTS KIT-MOTOR ASSM TRUTECH2097829-127 and associated hardware.	1, 2, 3, 4, 6	1, 3, 8	9
Drive DC Servo Kollmorgan PC2303-012-N 2097829-033 and associated hardware.	1, 2, 3, 4, 6	1, 3, 8	9
Drive DC Servo Quantum/Trutech PC2303-012-N 2097829-125 and associated hardware.	1, 2, 3, 4, 6	1, 3, 8	9
PART KIT-FIRMWARE 1.19.109 2097829-129	1, 2, 4	---	---
PART KIT-CONVERSION KIT TRUTECH/QUANTUM2097829-130	1, 2, 3, 4, 6	1, 3, 8	9
KIT T2100-ST CALIBRATION SOFTWARE 1_0_1_28 2097829-128	1, 2, 4	---	---
New Adjustment Hardware Kit 2097829-098	1, 2, 4	1, 6, 7, 8	9
Drive Board Assembly Hardware Kit 2097829-099	1, 2, 4	1, 6, 7, 8	9
Speed Sensor 3 Wire Hall Effect 2097829-075	1, 2, 3, 4	1, 2, 3, 8	9
Harness Speed Sensor 3 Wire 2097829-076	1, 2, 3, 4	1, 2, 3, 8	9
Tether clip with Magnet 2097829-120	1	2	---
BOARD FG0495-3 UNIVERSAL POWER SUPPLY FW V1.19.109 2097829-126	1, 2, 3, 4, 6	1, 3, 4, 8	9
Board FG0495-1 Universal Power Supply FW V1.66 2097829-022 and associated hardware.	1, 2, 3, 4, 6	1, 3, 4, 8	9
Board FG0495-3 Universal Power Supply FW V1.18.12 2097829-080 and associated hardware.	1, 2, 3, 4, 6	1, 3, 4, 8	9
Wheel Landing Gear 2097829-024	---	5	---
Pull Tether Assembly 2097829-026	1, 2, 3, 4	1, 2, 8	9
Emergency Stop with Guard Assembly 2097829-011	1, 2, 3, 4	1, 2, 8	9
Actuator 115V 2097829-028 or Actuator MSI 2097829-088 and associated hardware.	1, 2, 3, 4, 6	1, 4, 8	9
Actuator 220V 2097829-029 or Actuator MSI 2097829-089 and associated hardware.	1, 2, 3, 4, 6	1, 4, 8	9
Fan Assembly DC W/Connector 2097829-030 and associated hardware.	1, 2, 3, 4	1, 8	9

FRU Repairs / Exchange

Description	Tools	Functional Checkout Procedure	Electrical Safety Checks
Resistor Omite 85/150MJ40R0DE 209/829-031 and associated hardware.	1, 2, 3, 4	1, 8	9
Filter CorCom F4635/10EHT1 209/829-032/082 and associated hardware.	1, 2, 3, 4	1, 2, 8	9
Breaker 3 amp 209/829-034	1, 2, 3, 4	1, 2, 8	9
Breaker 1/4 amp 209/829-035	1, 2, 3, 4	1, 2, 8	9
Switch TA35-C324-L200C0 209/829-036	1, 2, 3, 4	1, 2, 8	9
Switch Mechanical Product 2421-305-119-2000 209/829-123	1, 2, 3, 4	1, 2, 8	9
Power Cord Inlet C-20 4793.4000 209/829-037	1, 2, 3, 4	1, 2, 8	9
Harness 110V DC Treadmill 209/829-040	1, 2, 3, 4	1, 2, 8	9
Harness 220V DC Treadmill 209/829-041	1, 2, 3, 4	1, 2, 8	9
Harness Resistor 209/829-042	1, 2, 3, 4	1, 2, 8	9
Speed Sensor Hall Effect 209/829-043	1, 2, 3, 4	1, 2, 3, 8	9
Harness Speed Sensor Hall Effect 209/829-044	1, 2, 3, 4	1, 2, 3, 8	9
Harness Elevation 209/829-045	1, 2, 3, 4	1, 2, 8	9
Harness Elevation FG0495-3 209/829-083	1, 2, 3, 4	1, 2, 8	9
Harness Motor Controller 209/829-046	1, 2, 3, 4	1, 2, 8	9
Harness Main DC Treadmill 209/829-047	1, 2, 3, 4	1, 2, 8	9
Harness Main DC Treadmill FG0495-3 209/829-084	1, 2, 3, 4	1, 2, 8	9
Relay 110V 209/829-085	1, 2, 3, 4	1, 2, 8	9
Relay 220V 209/829-086	1, 2, 3, 4	1, 2, 8	9
Surge Suppressor 209/829-048	1, 2, 3, 4	1, 2, 8	9
Cable RS232 Panel Mount 209/829-049	1, 2, 3, 4	1, 2, 8	9
Stem Bumper 209/829-052	1, 2	1, 2, 8	---
Bumper 1.8" Dia Black 209/829-053	1, 2	1, 2, 8	---
Belt Drive 209/829-054 and associated hardware.	1, 2, 3, 4	1, 3, 8	9
Belt Running with V-Guide 209/829-057 and associated hardware.	1, 2, 4	1, 6, 7, 8	9
Deck Isolator 1.75" Dia with Stud 209/829-058 and associated hardware.	1, 2, 4	1, 6, 7, 8	9
Deck Running Board 209/829-059 and associated hardware.	1, 2, 4	1, 6, 7, 8	9
Roller Front 209/829-060 and associated hardware.	1, 2, 4	1, 6, 7, 8	9

Functional / Post Repair / Preventative Maintenance Checkout Procedures

FRU Repairs / Exchange

Description	Tools	Functional Checkout Procedure	Electrical Safety Checks
Rear Roller 2097829-061 and associated hardware.	1, 2, 4	1, 6, 7, 8	9
End Cap LH 2097829-062 and associated hardware.	1, 2	1, 8	---
End Cap RH 2097829-063 and associated hardware.	1, 2	1, 8	---
Extrusion Anti-Skid Rubber 60" LG 2097829-064	1, 2, 4	1, 8	9
Handrail RH 2097829-066 and associated hardware.	1, 2, 3, 4	1, 2, 5, 8	9
Handrail LH 2097829-067 and associated hardware.	1, 2, 3, 4	1, 2, 5, 8	9
Handrail Center 2097829-068 and associated hardware.	1, 2, 3, 4	1, 2, 5, 8	9
Hood Kit 2097829-069	1, 2, 3, 4	1, 8	9
Test Plug 2097829-070	1, 2	1, 5, 8	---
SPARE COMM CABLE CLAMPS 2097829-074	1, 2	1	---
POWER CORD 2097829-101 through 2097829-107	1, 2, 4	1, 2, 8	9
USB CABLE 5M A TO B 2044095-001 (INTERFACE T2100-ST TO CARDIOSOFT/CS)	1, 2, 4	1, 8	---
CABLE RS232 INTERFACE T2100-ST TO MAC5500ST 2097829-108 (contains 2097829-008)	1, 2, 4	1, 8	---
CABLE RS232 INTERFACE T2100-ST TO CASE/PC 2097829-109 (contains 2097829-009)	1, 2, 4	1, 8	---
CABLE RS232 INTERFACE T2100-ST TO MAC2000ST 2097829-110 (contains 2097829-010)	1, 2, 4	1, 8	---
Final Checkout Preventative Maintenance Test			
Procedure (for PM performance)	Tools	Functional Checkout Procedure	Electrical Safety Checks
Preventative Maintenance Visual/Functional Checks			
Operation Check	1	2	---
Test Plug Check	1	5	---
Walking Belt Tension / Tracking	1, 2, 6	6	---
Walking Belt Speed Verification	1, 2, 6	7	---
Leakage Test	1, 2, 4, 5, 6	9	9

A

Serial # _____ Date Purchased _____

[illegible]

Reference Documents

If needed, the following installation documents can be referenced and may provide some helpful OS specific guidance for driver installation.

These files are PDF documents located in the directory called **FTDI PDF** on the flash drive shipped with the T2100-ST series treadmill.

- Windows 7:
AN_119_FTDI_Drivers_Installation_Guide_for_Windows_7.pdf
- Windows 8:
AN_234_FTDI_Drivers_Installation_Guide_for_Windows_8.pdf
- Windows 10:
AN_396_FTDI_Drivers_Installation_Guide_for_Windows_10.pdf

NOTE:

Do not download the driver from the FTDI website as referenced in Chapter 2 of these three documents

Reference Documents

Reprocessing Instructions

The following reprocessing instructions apply to all treadmills. The treadmill is neither sterile nor can it be sterilized.

Manufacturer: Full Vision Inc		Symbol: 
Device(s): Treadmills		
WARNINGS	Before cleaning the device, turn the main power switch to OFF and disconnect the treadmill from its power outlet. Never use wet cleaning materials near a power source: you could suffer an electrical shock. Do not improperly clean the device ○ During cleaning, do not allow fluid to enter the motor pan assembly area. ○ Always wear appropriate PPE when cleaning. ○ Do not use abrasive brushes or abrasive cleaners. This could damage the paint and plastic surfaces. ○ Do not use sharp tools (e.g. knife, metal scraper) or aggressive cleaning solvents for cleaning ○ Alcohols are flammable and should only be used in well ventilated spaces	
Limitations on Processing	Not Applicable	

Reference Documents


INSTRUCTIONS	
Initial Treatment at the Point of Use	Use a soft, clean, lint-free cloth / paper towel dampened (not wet) with tap water and mild detergent to wipe all treadmill surfaces to remove excess soil, moisture, and perspiration.
Preparation Before Cleaning	All cleaning solutions should be prepared at the dilution and temperature recommended by the manufacturer.
Cleaning: Automated	Not applicable
Cleaning: Manual	Step 1 – Turn the main power switch to OFF and disconnect the treadmill from its power outlet. Step 2 – Use 20 mL (0.68 oz) of an alkaline cleaner or similar mild, non-abrasive detergent (ex. Spray 409) on a soft, clean, lint-free 25 x 25 cm (9.8 x 9.8 in) size cloth / paper towel and manually clean item 1, 2. Reference Image 1 for corresponding items. Note: Do not spray cleaning chemical directly on the device Step 3 – Use an appropriately sized soft-bristled brush (e.g. standard toothbrush) dampened with 5 mL (0.17 oz) of mild soapy tap water (or similar mild, non-abrasive detergent) to clean item 3 and hard to reach spots. Reference Image 1 for corresponding items. Step 4 – Wipe all surfaces 1-5 (reference Image 1) with a soft, clean, lint-free 25 x 25 cm (9.8 x 9.8 in) size cloth / paper towel dampened with 20 mL (0.68 oz) (not wet) of lukewarm tap water 27°C to 44°C (80°F to 111°F) for a minimum of 30 seconds. Step 5 – Dispose of all cleaning materials according to your institutions established procedure.
Disinfection	If necessary, to disinfect the treadmill, follow these steps. Step 1 – Turn the main power switch to OFF and disconnect the treadmill from its power outlet. Step 2 – Prepare 20 mL (0.68 oz) of 70% or greater isopropyl alcohol disinfectant solution according to the manufacturer's directions. Step 3 – Use a soft, clean, lint-free 25 x 25 cm (9.8 x 9.8 in) size cloth / paper towel or gauze that is saturated with 20 mL (0.68 oz) of the disinfection solution and manually disinfect all contaminated surfaces of the device for a minimum of 15 minutes. Step 4 – Use an appropriately sized soft-bristled brush (e.g. standard toothbrush) dampened with 5 mL (0.17 oz) of the disinfectant solution to clean item 3 and hard to reach spots. (Reference Image 1 for corresponding items) Surface must remain wet for a minimum of 15 minutes. Step 5 – Wipe all contaminated surfaces 1-5 (reference Image 1) with a soft, clean, lint-free 25 x 25 cm (9.8 x 9.8 in) size cloth / paper towel or gauze dampened with 20 mL (0.68 oz) of purified water 27°C to 44°C (80°F to 111°F) for a minimum of 30 seconds. Step 6 – Dispose of all cleaning materials according to your institutions established procedure.

Reference Documents

Drying	Dry the device with a soft, clean, lint-free cloth/paper towel 25 x 25 cm (9.8 x 9.8 in).
Maintenance, Inspection and Testing	Visually examine each device for cleanliness. If visible soil remains, repeat cleaning procedure until the device is thoroughly clean.
Packaging	Not applicable
Sterilization	Not applicable
Storage	 The diagram shows three boxes for storage and operating conditions. The first box indicates '90% Non Condensing' humidity and '10%' relative humidity. The second box shows a temperature range from '-40°C' to '70°C' for 'Storage Temp.'. The third box shows a temperature range from '-5°C' to '38°C' for 'Operate Temp.'.
Additional Information	No additional requirements
Manufacturer Contact	Full Vision Inc. Email: tmservice@full-vision.com

Reference Documents

Image 1

Item	Description
1	Emergency Stop Button
2	Patient Grab Handrails
3	Pull Tether Switch
4	Hood
5	Running Belt

The instructions provided have been validated by the manufacturer of the medical device as being capable of preparing a medical device for reuse. It remains the responsibility of the processor to ensure that the processing, as actually performed using equipment, materials and personnel in the processing facility, achieves the DESIRED result. This requires verification and/or validation and routine monitoring of the process.

3017 Full Vision Drive
 Newton, KS 67114 USA
<http://www.full-vision.com>
 Phone: (316)-283-3344
 Fax: (316)-283-9522
 Email: sales2@full-vision.com

(Regulatory affairs only)
 Emergo Europe
 Westervoortsedijk 60
 6827 AT Arnhem
 The Netherlands

DISTRIBUTOR / Private Label Manufacturer:
 GE Medical Systems Information Technologies, Inc.

Distributor World Headquarters:

GE Medical Systems
 Information Technologies, Inc.
 9900 Innovation Drive
 Wauwatosa, WI 53226 USA

**GE Medical Systems Information Technologies, Inc., a General Electric Company,
 going to market as GE Healthcare.**

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
The screenshot shows a software window titled "system Configuration" with a close button in the top right corner. Below the title bar is a tabbed interface with tabs for "General", "Device", "Modem", "MUSE", "Option Code", and "Country Settings". The "Device" tab is selected. The window contains several configuration fields: "Treadmill Model" (dropdown menu showing "T2100"), "Port" (dropdown menu showing "COM 1"), "BP Monitor" (checkbox, unchecked), "BP Port" (dropdown menu showing "COM 1"), "Metabolic Cart" (checkbox, unchecked), "Engineer Model" (dropdown menu showing "internal"), "Port" (dropdown menu showing "COM A"), "SpO2 Monitor" (checkbox, unchecked), "SpO2 Port" (dropdown menu showing "COM 1"), "Ambulatory Blood Pressure" (checkbox, unchecked), "Ambulatory BP Device Port" (dropdown menu showing "COM 2"), "Sphygmomanometer Device Type" (checkbox, unchecked), "Sphygmomanometer Device Port" (dropdown menu showing "COM 2"), "Path to Holter ECG Program" (text field), "Path to External Program" (text field), "Use External Program for Sphygmomanometer" (checkbox, unchecked), "Procedure Name for External Program" (text field), "Analog 1" (dropdown menu showing "none"), "Analog 2" (dropdown menu showing "none"), "Analog 3" (dropdown menu showing "none"), "Analog 4" (dropdown menu showing "none"), "TTL output" (checkbox, unchecked), "Pulse Width (ms)" (text field showing "50"), "Delay (ms)" (text field showing "0"), and "Polarity" (dropdown menu showing "positive"). At the bottom of the window are buttons for "Print", "Help", "OK", and "Cancel".

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
The screenshot shows the "System Configuration" window with the "Devices" tab selected. The window has a title bar "System Configuration" and a close button. Below the title bar is a tab bar with "General", "Devices", "Modem", "MUSE", "Option Code", and "Country Settings". The "Devices" tab contains several settings: "Treadmill Model" (dropdown with "T2100" selected), "Port" (dropdown with "COM 1" selected), "BP Monitor" (dropdown with "no" selected), "Port" (dropdown with "COM 1" selected), "Metabolic Cart" (checkbox, unchecked), "Path to Hoker ECG Program" (text field), "Path to External Program" (text field), "Use External Program for Sponometry" (checkbox, unchecked), "Procedure Name for External Program" (text field), "Analog 1" (dropdown with "none" selected), "Analog 2" (dropdown with "none" selected), "Analog 3" (dropdown with "none" selected), "Analog 4" (dropdown with "none" selected), "Ergometer Model" (dropdown with "internal" selected), "Port" (dropdown with "COM A" selected), "SpO2 Monitor" (dropdown with "no" selected), "Port" (dropdown with "COM 1" selected), "Ambulatory Blood Pressure" (dropdown with "no" selected), "Ambulatory BP Device Port" (dropdown with "COM 2" selected), "Sponometry Device Type" (dropdown with "no" selected), "Sponometry Device Port" (dropdown with "COM 2" selected), "TTL output" (checkbox, unchecked), "Pulse Width (ms)" (text field with "50"), "Delay (ms)" (text field with "0"), and "Polarity" (dropdown with "positive" selected). At the bottom are buttons for "Print", "Help", "OK", and "Cancel".

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TANGO M2[®]

SERVICE MANUAL

SunTech Medical, Inc.
507 Airport Boulevard
Morrisville, North Carolina 27560-8200
United States of America
Toll free: 800.421.8626
Phone: 919.654.2300
Fax: 919.654.2301

SunTech Medical, Ltd.
Oakfields Industrial Estate
Stanton Harcourt Road
Eynsham, Oxon OX8 1TS
United Kingdom
Phone: (44)1865-884234
Fax: (44)1865-884235

Part # 80-0056-00 Rev A

MANUAL CONTENTS

SECTION 1	Introduction	3
SECTION 2	Inputs & Outputs	4
	Inputs	4
	Outputs	4
	Combined Inputs & Outputs	4
SECTION 3	Operation Modes & Menu System	5
	Normal Operation	5
	Diagnostics Menu	5
SECTION 4	Routine Maintenance	7
	Verification of Calibration	7
	Pressure Calibration	7
	Leak Rate Test	8
	Timeout Check	8
	Overpressure Check	9
SECTION 5	Non-Routine Maintenance	9
	General Internal Inspection	9
SECTION 6	Disassembly & Re-assembly Procedures	10
Appendix 1	Spare Parts List	12
Appendix 4	Key to Exploded Views	13
Schematics	PCB Schematics	CD
Drawings	Exploded View of Tango M2 Assembly	CD

This Manual and associated CD contains drawings and/or specifications that are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in strict confidence and shall not be reproduced or copied or used as the basis for manufacture or sale of apparatus without permission.

TANGO M2 STRESS TEST BLOOD PRESSURE MONITOR

1.0 INTRODUCTION

- 1.1 The SunTech Medical Tango M2 Stress Non-Invasive Blood Pressure Monitor is used to obtain blood pressure and heart rate readings during exercise testing. It can be used as a stand alone monitor or in conjunction with a host Stress Test system. The essential set-up consists of the monitor, pneumatic/K-sound patient cable and Orbit-K blood pressure cuff with microphone.**
- 1.2 The Tango M2 is typically interfaced to a host stress test system and receives its analog or digital ECG trigger via a coaxial BNC connection. Communication between the stress system and the monitor is via an RS-232 null modem cable or USB-A to USB-B cable. An optional internal ECG board is available for standalone operation.**

This Manual and associated CD contains drawings and/or specifications that are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in strict confidence and shall not be reproduced or copied or used as the basis for manufacture or sale of apparatus without permission.

2.0 TANGO M2 INPUTS & OUTPUTS

2.1 TANGO M2 INPUTS

- 2.1.1 **Power Switch.** located on rear of enclosure. This is a latching, push button on/off switch that removes or provides power to the monitor circuitry.
- 2.1.2 **External (Analog or Digital) ECG Trigger Input.** This is a BNC connection which provides the ECG signal from the host stress system to the Tango M2.
- 2.1.3 **Optional ECG Patient Cable (3 lead to 6-pin DIN).** Provides the ECG signal to the Tango M2. Used only if the Tango M2 is not using an external ECG source.
- 2.1.4 **Front Panel Keypad.** Six-key membrane pad used to control the Tango M2 operations.
- 2.1.5 **9VDC medical grade power supply, SunTech PN 19-0012-01.**

2.2 TANGO M2 OUTPUTS

- 2.2.1 **Cold cathode backlit LCD display.**
- 2.2.2 **Beeper (for audible feedback).**

2.3 TANGO M2 COMBINED INPUTS & OUTPUTS

- 2.3.1 **Pneumatic / K-Sound Patient Cable.** Connects the Tango M2 pneumatic pump to the arm cuff, and the Tango M2 K-Sound circuitry to the K-Sound microphone.
- 2.3.2 **9-Pin, Bi-directional RS-232 serial port.** Used for communicating with the host system.
- 2.3.3 **USB-B port used for communicating with the host stress system (alternative to the RS-232 serial port).**
- 2.3.4 **USB-A port used for downloading measurement history from the Tango M2 or uploading software/firmware updates.**

This Manual and associated CD contain data and information that are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in strict confidence and shall not be reproduced or copied or used as the basis for manufacture or sale of apparatus without permission.

3.0 OPERATIONAL MODES & MENU SYSTEM DESCRIPTION

- 3.1 **Normal Operation** - The Tango M2 is a menu-driven system. All settings and operation are controlled by using the front panel keypad. Refer to the Tango M2 User's Manual for detailed instructions on using the system under normal operational conditions.
- 3.2 **Diagnostics Menu** - Certain maintenance operations may be accessed by a hidden key sequence that is not normally available to the user. The Maintenance Menu is accessed by the following keypad input sequence:

- 1. Press and hold the MENU key.
- 2. Press and release the UP arrow key.
- 3. Press and release the Down arrow key.
- 4. Press and release the Down arrow key.
- 3. Press and release the UP arrow key.
- 4. Release the Menu key.

The available options are shown on-screen and are more fully discussed in the following sections.

In addition to the options used by the standard maintenance procedures (described later in this manual), the Diagnostics Menu also provides access to the Configuration Options settings and Direct Hardware Control.

- 3.2.1 **Verify Calibration** - Selecting this option will display the current pressure within the Tango M2 system. This option should be used with a T-Tube kit (SunTech part number 98-0030-00). Overpressure and safety timeout (pressure above 15mmHg in system for greater than 180 seconds) are disabled. Also, the system will display values up to 320mmHg. Exiting the diagnostic menu will enable overpressure, the safety timeout and the display will no longer display pressures above specified ranges as provided in the User Manual.

3.2.2 Calibrate

- 3.2.2.1 **Verify Calibration** - Selecting this option will display the current pressure within the Tango M2 system. This option should be used with a T-Tube kit (SunTech part number 98-0030-00). Overpressure and safety timeout are enabled. The monitor will only display pressures up to maximum specified ranges as provided in the User Manual.

- 3.2.2.2 **Calibrate** - **SHOULD ONLY BE PERFORMED BY QUALIFIED PERSONNEL USING A CALIBRATED REFERENCE PRESSURE! THIS SET-UP IS CRITICAL TO THE ACCURACY OF THE TANGO M2 MONITOR!** The person calibrating the monitor must allow the pressure to stabilize at each pressure point setting before pressing ENTER. Two pressure points will be set during the calibration, 0mmHg and 250mmHg. When complete, Verify

This Manual and associated CD contains drawings and/or specifications that are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in strict confidence and shall not be reproduced or copied or used as the basis for manufacture or sale of apparatus without permission.

Calibration must be performed. See Section 4.3 for more details on performing the calibration.

3.2.3 Test Input Channels – SunTech use only.

3.2.4 Test Outputs

3.2.4.1 Pump – Pumps run at full power when set to On and shut off fully when set to OFF.

3.2.4.2 Valve – Drop valve closes when set to Closed and opens when set to Open.

3.2.4.3 Valve – Step valve closes when set to Closed and opens when set to Open.

3.2.4.4 Control Valve – Linear control valve closes completely when set to Closed and opens completely when set to Open.

3.2.4.5 Beeper – sounds Beeper when set to On and turns off beeper when set to OFF.

3.2.5 Configuration

3.2.5.1 ECG Skip – This setting is used to help determine the K-sound window from ECG. This value is normally set to 3 when the ECG trigger is set to Digital. The value is set to 0 when the ECG trigger is set to Analog. Controls are Auto (0 or 3, depending on ECG trigger setting) and 1 – 15 in increments of 1 and will override the aforementioned values associated with the ECG trigger setting. **Note: Any change in this value can cause incorrect detection of K-sounds.**

3.2.5.2 Delay – This delays the window of K-sound detection and should normally be set to 0ms. This setting is used to compensate for any delay that the host stress test system might induce in the ECG trigger. Delays range from Default (no delay) to 600ms in increments of 10ms. **Note: Any change in this value can cause incorrect detection of K-sounds.**

3.2.6 Transition Mode – SunTech use only.

3.2.7 Special

3.2.7.1 Dump Data to USB – In the event of a system error by the Tango M2, this utility will dump debug information to a removable Flash Drive (USB 2.0) for analysis by SunTech.

3.2.7.2 Service Menu – SunTech use only.

4.0 TANGO M2 ROUTINE MAINTENANCE PROCEDURES

- 4.1 The following standard maintenance procedures are provided to ensure correct calibration of the Tango M2 and correct operation of the basic patient safety systems of the Tango M2. SunTech recommends that the calibration and safety systems be checked annually.

4.2 VERIFICATION OF PRESSURE CALIBRATION

Equipment Required:

1. Calibrated reference sphygmomanometer
2. 500ml volume.
3. Hand Inflation Bulb with bleed valve
4. Tubing, tee pieces and miscellaneous connectors.

- 4.2.1 Tee the sphygmomanometer, the 500ml volume and the inflation bulb together and connect to the patient cable hose connection.
- 4.2.2 Within the Diagnostics Menu, select Verify Calibration. The Tango M2 will now have its valves closed and will display the pressure applied to the patient hose connector.
- 4.2.3 Verify the Tango M2 calibration by inflating manually and checking the sphygmomanometer against the pressure reading on the Tango M2 display. Verify the readings as shown in Table 1 below.

Calibrated Reference Pressure (mm Hg)	Tango M2 Displayed Pressure (mm Hg)
0	0
50	48 - 52
100	98 - 102
150	148 - 152
200	198 - 202
250	248 - 252
280 (note)	278 - 282

Table 1: Pressure Verification Table

NOTE – Increase pressure to 280mmHg slowly, to avoid overshoot tripping the overpressure limit.

- 4.2.4 If any pressure reported by the Tango M2 is outside limits, the monitor must be calibrated (see section 4.3).
- 4.2.5 Press ENTER to open the valves and release the pressure from the system.
- 4.2.6 Selecting ENTER or pressing the Start/Stop button will leave the Diagnostic Menu.

4.3 SET CALIBRATION

This Manual and associated CD contains drawings and/or specifications that are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in strict confidence and shall not be reproduced, copied or used as the basis for manufacture or sale of apparatus without permission.

Equipment Requirements:

1. Calibrated aneroid sphygmomanometer
2. 500ml volume
3. Hand inflation bulb with check valve
4. Tubing, syringe pump and standard luer lock connectors

- 4.3.1 Tee the sphygmomanometer, the 500ml volume and the inflation bulb together and connect to the patient cable hose connection.
- 4.3.2 Access the **Flag** in the **Flags** Menu (see section 3.2).
- 4.3.3 Select **Calibrate** in the **Diagnostics Menu** and select **Calibrate** again.
- 4.3.4 Operate the hand bulb on the pneumatic bulb to obtain 0mmHg pressure. Ensure the aneroid sphygmomanometer displays 0mmHg and then press **ENTER** to store the calibration in the Tango M2.
- 4.3.5 Close the check valve on the bulb and inflate to 250mmHg. Allow the pressure to stabilize as necessary to maintain the pressure at exactly 250mmHg on the aneroid sphygmomanometer. Press **ENTER** to set the 250mmHg point on the Tango M2. **Note:** A failure will be displayed if the monitor recognizes a pressure outside reasonable expected limits of a clinical use. If a failure is indicated, the user will be prompted to cancel (that can be done by pressing the **ESC** key).
- 4.3.6 Check the pressure calibration.
- 4.3.7 Select **Verify** in the **Diagnostics Menu** steps 4.2.3 through 4.2.6. In the event verification fails, repeat every step in section 4.3. If verify calibration continues to fail, do not use the Tango M2 monitor and contact SunTech Medical.

4.4 LEAK TEST PROCEDURE

Equipment Requirements:

1. Stop Valve
2. 500ml volume
3. Hand inflation bulb with check valve
4. Tubing, syringe pump and standard luer lock connectors.

- 4.4.1 Tee the inflation bulb, the 500ml volume and the inflation bulb together and connect to the patient cable hose connection.
- 4.4.2 Select **Calibrate** in the **Diagnostics Menu**. The Tango M2 will now be set to the 0mmHg point. Inflate the pressure applied to the patient hose connection.
- 4.4.3 Inflate the pressure to 250mmHg. Allow the system to stabilize for 1 minute, stop the inflation and allow the pressure to fall. Note the exact pressure, and allow the pressure to stabilize again. After 60 seconds have elapsed, the pressure should not be any greater than 6mmHg.
- 4.4.4 If the pressure is greater than 6mmHg, increase pressure from the system.

4.5 SAFETY CHECK PROCEDURE

Equipment Requirements:

1. Stop Valve
2. 500ml volume

This Manual and associated CD-ROM are the property of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. © 2004 SunTech Medical, Inc. All rights reserved. No part of this manual may be reproduced or copied or used as the basis for manufacture or distribution without the written permission of SunTech Medical, Inc.

3. Hand Inflation Bulb with Head Valve
4. Tubing, Tee pieces and miscellaneou s
- 4.5.1 Tee the 500ml volume and the inflation bulb together and connect to the patient cable (see Figure 4.5.1).
- 4.5.2 Select Calibrate from the menu and then select Verify Calibration. Note: It is important to select Verify Calibration from the Calibrate menu as this has to be enabled.
- 4.5.3 Using the inflation bulb, inflate the system to approximately 200mmHg per the monitor's display. Start the stopwatch as soon as pressure exceeds 180mmHg. Release the pressure from the system after 15 to 18 seconds. Stop the stopwatch as soon as pressure falls below 15mmHg. This time must be 10 seconds or less.

4.6 OVERPRESSURE CHECK

Equipment Required:-

1. Stop Watch.
2. 500ml volume.
3. Hand Inflation Bulb with Head Valve
4. Tubing, Tee pieces and miscellaneou s
- 4.6.1 Tee the 500ml volume and the inflation bulb together and connect to the patient cable (see Figure 4.6.1).
- 4.6.2 Select Calibrate from the menu and then select Verify Calibration. Note: It is important to select Verify Calibration from the Calibrate menu as this has to be enabled.
- 4.6.3 Using the inflation bulb, inflate the system to approximately 280mmHg pressure per the monitor's display. Start the stopwatch to inflate the system until the UUT opens the bulb to release the pressure. This must occur between 290mmHg and 300mmHg.

5.0 TANGO M2 NON-ROUTINE MAINTENANCE PROCEDURES

5.1 Internal Maintenance

WARNING – There are no user serviceable parts within the Tango M2. However, the enclosure may be removed and the components separated if it is required to inspect the unit for mechanical damage due to foreign ingress after spillage, or if the unit has been dropped.

5.2 General Internal Inspection

WARNING – The warranty for the Tango M2 will be voided if the warranty seal is broken. Contact SunTech Medical for service is required during the warranty period. The Tango M2 enclosure should only be opened by qualified personnel.


- 5.2.1 Open the enclosure by removing the four Torx head screws located on the four corners of the rear enclosure (one corner has the warranty seal).
- 5.2.2 See section 6 for disassembly instructions.
- 5.2.3 Inspect assembly for damage or contamination. Should any damage or liquid contamination be evident, the Tango M2 should be re-assembled and returned to SunTech for board replacement. Any dry contamination

This Manual and associated CD contains drawings and/or specifications of SunTech Medical, Inc; 507 Airport Blvd; Suite 117; Morrisville, NC 27560; USA. Copies are issued in the English language. This manual is not to be used as the basis for manufacture or repair of equipment.

5.2.5 Perform a Leak Test and a Timeout Test (as described in Section 4) to verify that the purifier system is functioning correctly. If either fail, open the monitor back up and check for damage and kinks in the pneumatic tubing.

K-Sound Cable

Pneumatic Tubing

Keypad Cable

- 6.1 Remove the four screws from the front enclosure. The gasket seal must be broken to avoid the warranty on the front enclosure.
- 6.2 Remove the front enclosure away from the device.
- 6.3 Remove the UI module cable from the UI module. Lay the UI module facing up. Remove the UI module assembly to the rear.
- 6.4 Remove the UI module of the Tango M2 board.
- 6.5 Remove the manifold and K-sound completely removed from the rear.

- 6.5.1 Pneumatic manifold: Connect the pneumatic line connection on the rear enclosure. Thread the pneumatic line through the I/O board so that the tubing exits straight out of the back of the unit.
- 6.5.2 K-sound cable: Disconnect the K-sound cable from J1 on the DKA module located underneath the unit. Plug the cable back through the I/O board. The connector on the cable should be inserted into the rear enclosure.

REASSEMBLY

- 6.6 Place the rear enclosure on the unit. The pneumatic manifold, K-sound cable, and power button must be connected first before the unit can be fully inserted back into the enclosure.
 - 6.6.1 K-sound cable: Thread the K-sound cable through the larger hole in the I/O board and connect the cable to J1 on the DKA module.
 - 6.6.2 Pneumatic manifold: Connect the pneumatic line through the rectangular cutouts on the I/O and UI boards. Plug the pneumatic line into the barbed connector.
 - 6.6.3 Power button: Slide the power button into the power switch.
 - 6.6.4 Place the rear enclosure on the unit. Insert the unit onto the four guide pins. Ensure that the cable is not crimped by the K-sound cable or by the unit itself. The cable should not form a loop.
- 6.7 Secure the module to the rear enclosure.
 - 6.7.1 Insert a screw and washer into the rear enclosure. Tighten securely.
- 6.8 Place the Targeted unit on the rear enclosure. Connect the front enclosure keypad cable to the UI module.
- 6.9 Place the Targeted ML display on the rear enclosure. Insert the four screws located in the slots on the back of the display. Tighten securely, ensuring that the keypad cable is not crimped.
- 6.10 Perform a Leak Rate test and a Functional Test (as described in Section 4) to verify that the pneumatic system is functioning properly. If the test fails, open the monitor back up and check for damage to the pneumatic system.

APPENDIX

APPENDIX 1

PART NUMBER	DESCRIPTION	QUANTITY	UNIT	REMARKS
91-0129-00	...	...	...	...
36-0001-01	...	...	...	...
98-0062-25	...	...	...	...
98-0235-00	...	...	...	...
98-0235-01	...	...	...	...
92-0138-00*	...	...	...	...
92-0139-00*	...	...	...	...
92-0140-00	...	...	...	...
98-1009-00	...	...	...	...
91-0004-00	...	...	...	...
19-0012-01	...	...	...	...
58-0029-01	...	...	...	...
55-0119-00	...	...	...	...
55-0095-00	...	...	...	...
55-0043-05	...	...	...	...
55-0143-00	...	...	...	...
17-0020-00	...	...	...	...
26-0012-00	...	...	...	...
39-0112-00	...	...	...	...
39-0101-00	...	...	...	...
39-0102-00	...	...	...	...
39-0103-00	...	...	...	...
12-0002-00	...	...	...	...

* Not available

This manual and associated software are the property of the U.S. Government and are not to be distributed outside the U.S. Government.

This manual and associated software are the property of the U.S. Government and are not to be distributed outside the U.S. Government.

This manual and associated software are the property of the U.S. Government and are not to be distributed outside the U.S. Government.

APPENDIX 2

Drawing Reference (97-0144-XX-A)	Quantity	Part Number
FH1	1	39-0102-00
LP1	1	39-0101-00
LENS1	1	58-0029-01
DISP1	1	26-0012-00
PWA2	1	92-0139-00
PWA1	1	92-0138-00
PB1	1	39-0102-00
M1	1	55-0119-00
MT1	1	39-0094-00
LBL1	1	71-0369-00
FT1	1	42-0001-00
RH1	1	39-0106-00

This Manual and associated CD contain
Morrisville, NC 27560, U.S.A. Copyright


117
for manufacture or

Note: - On production boards, footprints for the Micro and Standard size SD Card components will be etched but only one (or neither) of these will be populated

SD Card #	Codec	REV	CN NUMBER	DATE	APPROVAL
		A	4302	20 Feb 2012	K. Bigda
		B	4484	02 Aug 2012	K. Bigda


J1-A
300 SODIMM



J1-C
PX300 SODIMM

J2
n/f
uSD Card Socket

C2
n/f


To Codec


Suntech Medical Inc.
507 Airport Boulevard, Morrisville, NC 27560 (919) 654 2300


Title Tango MD UI Board

Drawing Number	92-0139-XX-SD	Rev	B
Creation Date	16 March 2011	Last Update	Thu Aug 02, 2012
File	Drawn by	A. Leech	Sheet 1 of 9


This drawing or specification is the property of SunTech Medical Inc. Copies are made for internal use only and shall not be reproduced, stored, or used as the basis for any other drawing or specification without express written permission by SunTech Medical Inc.



2017 Airport, vol. 8117, Appendix No. 2, 560 (919) 654-4 (10141)			
TANG DM2 MONITOR ASSEMBLY			
WKS. NO.	97-0144	REV	A
1	3	1	SWET 2 CM 2



 **ООО ЦПМИ**  **INFO@CPMI.RU**  **WWW.GOSZDRAVNADZOR.RU**

J1-I
PX300 SODIMM

SDA1 O-56

SCL1 O-53

SDA2 O-66

SCL2 O-69

J1-J
PX300 SODIMM

GPIO4 O-97

GPIO22 O-100

GPIO40 O-98

GPIO41 O-101

GPIO42 O-103

RS-223

RS-223


RS-223

Sup ch Mc

Rev 1.0

Ver 1.0


Sheet 9 of 9


Smith Medical Inc.
 17400 Rte. 100, Raleigh, NC 27606 (919) 854-2300
 1992 Model 50000, 50001, 50002, 50003, 50004, 50005, 50006, 50007, 50008, 50009, 50010, 50011, 50012, 50013, 50014, 50015, 50016, 50017, 50018, 50019, 50020, 50021, 50022, 50023, 50024, 50025, 50026, 50027, 50028, 50029, 50030, 50031, 50032, 50033, 50034, 50035, 50036, 50037, 50038, 50039, 50040, 50041, 50042, 50043, 50044, 50045, 50046, 50047, 50048, 50049, 50050, 50051, 50052, 50053, 50054, 50055, 50056, 50057, 50058, 50059, 50060, 50061, 50062, 50063, 50064, 50065, 50066, 50067, 50068, 50069, 50070, 50071, 50072, 50073, 50074, 50075, 50076, 50077, 50078, 50079, 50080, 50081, 50082, 50083, 50084, 50085, 50086, 50087, 50088, 50089, 50090, 50091, 50092, 50093, 50094, 50095, 50096, 50097, 50098, 50099, 50100

EC Channels


Optional Internal ECG Circuit - At Risk

16 Bit ADC

Medical Inc.
201 654 2300
Rev A Sheet 4 of 11

12	12	12
12	12	12
73PW	73PW	73PW
CC	CC	CC
08	08	08
01	01	01
02	02	02
03	03	03
04	04	04
05	05	05
06	06	06
07	07	07
LE	LE	LE

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
С. 1000-1100 Sheet 6 of 11



H051 111 Circuit

R34
100K

R34
100K

± K Flags
from micro

Medical In:
K5004
K5005
K5006
GND
VCC

Safety Micro

Reset
Circuit

1. 10k
2. 10k
3. 10k
4. 10k
5. 10k
6. 10k
7. 10k
8. 10k
9. 10k
10. 10k
11. 10k
12. 10k
13. 10k
14. 10k
15. 10k
16. 10k
17. 10k
18. 10k
19. 10k
20. 10k
21. 10k
22. 10k
23. 10k
24. 10k
25. 10k
26. 10k
27. 10k
28. 10k
29. 10k
30. 10k
31. 10k
32. 10k
33. 10k
34. 10k
35. 10k
36. 10k
37. 10k
38. 10k
39. 10k
40. 10k
41. 10k
42. 10k
43. 10k
44. 10k
45. 10k
46. 10k
47. 10k
48. 10k
49. 10k
50. 10k
51. 10k
52. 10k
53. 10k
54. 10k
55. 10k
56. 10k
57. 10k
58. 10k
59. 10k
60. 10k
61. 10k
62. 10k
63. 10k
64. 10k
65. 10k
66. 10k
67. 10k
68. 10k
69. 10k
70. 10k
71. 10k
72. 10k
73. 10k
74. 10k
75. 10k
76. 10k
77. 10k
78. 10k
79. 10k
80. 10k
81. 10k
82. 10k
83. 10k
84. 10k
85. 10k
86. 10k
87. 10k
88. 10k
89. 10k
90. 10k
91. 10k
92. 10k
93. 10k
94. 10k
95. 10k
96. 10k
97. 10k
98. 10k
99. 10k
100. 10k

10k

10k

Power Control

1. 10k
2. 10k
3. 10k
4. 10k
5. 10k
6. 10k
7. 10k
8. 10k
9. 10k
10. 10k
11. 10k
12. 10k
13. 10k
14. 10k
15. 10k
16. 10k
17. 10k
18. 10k
19. 10k
20. 10k
21. 10k
22. 10k
23. 10k
24. 10k
25. 10k
26. 10k
27. 10k
28. 10k
29. 10k
30. 10k
31. 10k
32. 10k
33. 10k
34. 10k
35. 10k
36. 10k
37. 10k
38. 10k
39. 10k
40. 10k
41. 10k
42. 10k
43. 10k
44. 10k
45. 10k
46. 10k
47. 10k
48. 10k
49. 10k
50. 10k
51. 10k
52. 10k
53. 10k
54. 10k
55. 10k
56. 10k
57. 10k
58. 10k
59. 10k
60. 10k
61. 10k
62. 10k
63. 10k
64. 10k
65. 10k
66. 10k
67. 10k
68. 10k
69. 10k
70. 10k
71. 10k
72. 10k
73. 10k
74. 10k
75. 10k
76. 10k
77. 10k
78. 10k
79. 10k
80. 10k
81. 10k
82. 10k
83. 10k
84. 10k
85. 10k
86. 10k
87. 10k
88. 10k
89. 10k
90. 10k
91. 10k
92. 10k
93. 10k
94. 10k
95. 10k
96. 10k
97. 10k
98. 10k
99. 10k
100. 10k

10k

10k

10k

10k

10k

10k

10k

10k

10k

10k

10k

Figure 1

Key-see
41943-1

47.15.5

**BP Micro - Pump & Valve
Gate Drive signals**

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DB
FS10

5

 $\leq$


15

200.2

2000

Medical In-

$$\frac{q_2 q_3 \dots q_n}{q_2 q_3 \dots q_n} = \frac{1}{1} = 1$$

Suntech Medical Inc.
507 Airport Boulevard, Morrisville, NC 27560 (919) 654 2300
Drawing Number 82-0133-XX-SD Rev A Sheet 11 of 11

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All Rights Reserved