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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of GUAN YINI of SHENZHEN WELL ELECTRIC COMPANY LIMITED and the seal of the said company on the annexed DOCUMENT are genuine.

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Authorized
Signature:

Zhou Siyuan

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(Date: Oct. 17, 2023)

证明书 址 Website for verification certificate: <http://www.rzccpit.com/validate.html>

Breath Alcohol Analyzers
АЛКОТЕКТОР PRO-100 touch-M

Operation Manual
РЭ 1.2

Sheets 88

APPROVED

Director

Shenzhen Well Electric Company Limited

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
2023

CONTENT

1 DESCRIPTION AND OPERATION	6
1.1 General provisions	6
1.1.1 Purpose	6
1.1.2 Scope	6
1.1.3 Designated users	6
1.1.4 Contraindications for use	6
1.1.5 Residual risk assessment	6
1.1.6 Operating conditions	6
1.1.15 Type of contact with the human body	7
1.2 Metrological and technical characteristics	7
1.2.16 Electrical power supply for analyzer	10
1.2.17 Electrical power supply of thermal printer	11
1.3 Composition of the analyzer	12
1.4 Design and operation of the analyzer	15
1.5 Marking and sealing	22
1.6 Packaging	23
2 INTENDED USE	23
2.1 Operating limitations	23
2.2 Turning on the analyzer	25
2.3 Analyzer menu	26
2.4 Functions of virtual touch screen buttons	28
2.5 Settings menu	29
2.5.5 Setting the measurement mode	30
2.5.6 Auto power off settings	31
2.5.7 Print Settings	32
2.5.8 Setting up Bluetooth connection	33
2.5.9 Adjusting the date and time	34
2.6 Preparing the analyzer for operation	35
2.6.1 Putting the analyzer into operation	35
2.6.2 Preparing the analyzer before starting work	36
2.7 Preparation and operation of the thermal printer	37
2.8 Operating procedure in measurement mode with data entry	39
2.9 Manual sampling mode	47
2.10 Refusal to measure	48
2.11 Operating procedure in measurement mode without data entry	49
2.12 Operating procedure in screening mode	50

2.13 Viewing and printing (transferring to an external device) measurement protocols saved in the analyzer's memory	53
2.14 Contents of the printed measurement report	56
2.15 Turning off the analyzer	57
2.16 Turning off the thermal printer	57
2.17 Possible malfunctions of the analyzer and thermal printer	58
3 MAINTENANCE	61
3.4 On-site maintenance	62
3.4.2 Charge/replace analyzer batteries	62
3.4.3 Charging the battery pack (for thermal printer)	63
3.4.4 Refilling the thermal printer with thermal paper	63
3.4.5 Touch screen calibration	64
3.4.6 Cleaning and disinfection	65
3.5 Periodic maintenance	65
4 TRANSPORTATION, STORAGE AND DISPOSAL	67
4.1 Transportation	67
4.2 Storage	67
4.3 Disposal	68
5 CHECKING THE ANALYZER	69
6 MANUFACTURER'S WARRANTIES AND REPAIRS	70
6.1 Manufacturer's warranty	70
6.2 Warranty repairs	70
6.3 Post-warranty maintenance and repair	70
APPENDIX A Regulatory documents	71
APPENDIX B Information on electromagnetic compatibility	72
APPENDIX C Weight and dimensions	76
APPENDIX D Signs and symbols used in marking	78
APPENDIX E Maintenance Tools	80
CHANGE REGISTRATION SHEET	82

This Operating Manual (hereinafter referred to as the OM) is intended to study the design, principle of operation, metrological and technical characteristics of the Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M (hereinafter referred to as the analyzer) and contains information necessary for its proper operation and technical service.

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Address : 1-3F, Building 101, No. 407, HeDongCun , HengKeng , GuanCheng Community, GuanHu Street, Longhua District, Shenzhen City, Guangdong, China, Tel: 86-755-83160728, Fax: 86-755-83160467, wellzp@well-co.com.

Authorized representative of the manufacturer (manufacturer):
Limited Liability Company "AJKOTEKTOP" (LLC "ALKOTEKTOR").

Location address: 191036, St. Petersburg, st. 1st Sovetskaya, 10, letter A, room . 2-N, tel. (812) 320-22-97 , email . address: www.alcotector.ru , email mail: info@alcotector.ru .

AJKOTEKTOP[®] is a registered trademark of ALKOTEKTOR LLC.

Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M is registered by the Federal Service for Surveillance in the Sphere of Health and is approved for circulation in the Russian Federation, registration certificate No. _____ dated _____.

Breath Alcohol Analyzers AJKOTEKTOP PRO-100 touch-M are registered with the Federal Information Foundation for Ensuring the Uniformity of Measurements under No. 56632-14.

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- ambient temperature, ° C: from 0 to plus 40;
- relative humidity of the ambient air (without moisture condensation), %: from 10 to 100;
- atmospheric pressure, kPa: from 84.0 to 106.7.

 *Note – To ensure safety when operating the analyzer with a power adapter connected, the requirements for the operating ambient temperature range of the power adapter must also be taken into account.*

1.1.7 Abbreviated designation of the analyzer: ALKOTEKTOR PRO-100 touch - M.

Note – On the front panel of the analyzer there is an abbreviated designation of the medical device in the form of an inscription "ALKOTEKTOR® PRO-100 touch - M".

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1.2.1 The measurement range and permissible error limits of the analyzer at ambient temperatures from plus 15 °C to plus 25 °C inclusive are given in Table 1.

Table 1 – Metrological characteristics of the analyzer at temperatures from plus 15 °C to plus 25 °C inclusive

Measuring range of ethanol mass concentration, mg/l	Limits of permissible error at temperatures from +15 to +25 °C incl.	
	absolute	relative
from 0 to 0.50 incl.	±0.05 mg/l	–
St. 0.50 to 0.95	–	±10%

Note – The analyzer has a programmatically set minimum interval of readings that are displayed on the analyzer screen and on paper in the form of zero readings: from 0.00 to 0.05 mg/l.

1.2.2 Indication range, mg/l: from 0.00 to 3.00.

Notes:

1 When the analyzer readings exceed the upper measurement limit, the mass concentration of ethanol in the selected sample of exhaled air exceeds 0.76 mg/l.

2 When the analyzer readings exceed the upper reading limit, the message ">3.00 mg/l" is displayed on the analyzer screen and in the printed measurement protocol (or transmitted to an external device).

1.2.3 Price of the lowest digit of the scale, mg/l: 0.01.

1.2.4 The limits of permissible error of the analyzer depending on the ambient temperature are given in Table 2.

Table 2 – Metrological characteristics of the analyzer depending on the ambient temperature

Ambient temperature	Measuring range of ethanol mass concentration, mg/l	Limits of permissible errors ¹⁾	
		absolute	relative
from 0 °C to +5.0 °C incl.	from 0 to 0.25 incl.	±0.05 mg/l	–
	St. 0.25 to 0.95	–	±20 %
St. +5.0 °C to +15.0 °C incl.	from 0 to 0.33 incl.	±0.05 mg/l	–
	St. 0.33 to 0.95	–	±15 %
St. +15.0 °C to +25.0 °C incl.	from 0 to 0.50 incl.	±0.05 mg/l ²⁾	–
	St. 0.50 to 0.95	–	±10% ²⁾
St. +25.0 °C to +40.0 °C incl.	from 0 to 0.50 incl.	±0.05 mg/l	–
	St. 0.50 to 0.95	–	±10 %

¹⁾ The table indicates the limits of permissible error of the analyzer under operating conditions given in 1.1.6 of this manual.
²⁾ According to table 1.

1.2.5 additional errors from the presence of unmeasured components.

1.2.6 Parameters of the analyzed gas mixture when a sample is supplied to the analyzer input (automatic sampling mode):

- flow rate of the analyzed gas mixture, l/ min: not less than 20;
- sample volume of the analyzed gas mixture, l: not less than 1.2.

1.2.7 Preparation time for work after inclusions, s: not more than 5.

1.2.8 Measurement time after sampling, s: not more than 15.

1.2.9 Time to prepare for work after measurement, s: not more than 20.

1.2.10 The analyzer provides ²the ability to print measurement reports on a thermal printer in a specified number of copies or transfer the measurement report to external devices at a distance of up to 10 m.

Bluetooth transmitter characteristics :

- frequency, MHz: from 2400 to 2483.5;
- modulation method: GFSK, 0.5 Gaussian ;
- radiation power, dBm , no more than: 4.

1.2.11 Printing of measurement reports is carried out on thermal paper installed in a thermal printer. The overall dimensions of the thermal paper roll are given in Appendix B.

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– mode for adjusting and checking readings;

– dry gas test mode – mode for measuring the mass concentration of ethanol vapor, which is used only to check the readings and calibrate the analyzer using ethanol/nitrogen gas mixtures in pressure cylinders;

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1.2.14 The analyzer's memory stores the results of the last 10,000 measurements.

1.2.15 Automatic shutdown time of the analyzer (set by the user), min: 1/2/5/10.

1.2.16 Electrical power supply for the analyzer

1.2.16.1 The analyzer can be powered electrically from:

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– power adapter for the analyzer with the characteristics specified in Table 3, connected to the AC mains;

– from the analyzer's on-board power adapter with the characteristics specified in Table 3, connected to the vehicle's on-board network;

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Analyzer on-board power adapter	Output voltage	≈ 5 V
---------------------------------	----------------	---------------

1.2.16.2 Number of measurements on the analyzer without replacing/charging batteries:

- without replacing batteries – at least 3000;
- without charging batteries – ^{at} least 1500 (when using batteries with a capacity of at least 2200 mAh).

1.2.17 Electrical power supply for thermal printer

1.2.17.1 The thermal printer can be powered electrically from:

- a battery pack with a voltage of 7.4 V and a capacity of at least 1800 mA · h;
- a charger power adapter ⁴ with the characteristics specified in Table 4, connected to an alternating current network;
- from the on-board power adapter of the charger ⁴ with the characteristics indicated in Table 4, connected to the vehicle's on-board network.

Table 4 - Characteristics of charger power adapters

Adapter type	Charger Power Adapter Specifications	
	parameter	meaning
AC adapter charger	Rated input voltage	~ 100-240 V
	Frequency	50/60 Hz
	Output voltage	≈ 12 V
	Output current	3 A
	Electric shock protection class	II
	Operating ambient temperature range	0°C to 40°C
On-board charger power adapter	Input voltage	≈ 12 V ⁵
	Output voltage	≈ 12 V

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⁵ The input voltage is determined by the model of the on-board adapter included in the delivery package of the analyzer; it can also be ≈ 24 V, $\approx 12-24$ V (indicated on the marking of the on-board power adapter).

1.2.17.2 The battery pack can be charged:

- in the power compartment of a thermal printer ⁶ when connecting a network or on-board power adapter to the charger;
- in a charger with the characteristics specified in Table 5, when connecting to it a network or on-board power adapter of the charger with the characteristics specified in Table 4.

Table 5 - Charger characteristics

Characteristic name	Meaning
Input voltage	≈ 12 V
Output voltage	≈ 8.4 V
Output current, no more	1 A

1.2.18 Overall dimensions of the analyzer (LxWxH), mm, no more: 150x70x35.

1.2.19 Weight of the analyzer (with batteries , without thermal printer), g, no more than: 280.

1.2.20 Weight and dimensional characteristics of components and accessories are given in Appendix B.

1.2.21 Service life of the electrochemical sensor installed in the analyzer, years: at least 2.

1.2.22 Operating time interval of the analyzer without adjusting the readings, months: not less than 12.

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1.2.23 Average service life of the analyzer, years: 5.

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2. Mouthpiece according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "ALCOTEKTOR mouthpiece with two holes" – no more than 5000 pcs.

3. USB cable – 1 pc.

4. Thermal printer – 1 pc.

5. Battery pack – no more than 2 pcs.

6. Thermal paper – no more than 20 pcs.

7. Charger – 1 pc.

8. Charger power adapter – 1 pc.

9. On-board charger power adapter – 1 pc.

10. Analyzer power adapter – 1 pc.

11. Analyzer on-board power adapter – 1 pc.

12. Analyzer power cable – 1 pc.

13. Interface cable – 1 pc.

1.3.2 The analyzer delivery set is shown in Table 6.

Table 6 – Analyzer contents

Name	Designation	Quantity
Ethanol vapor analyzer in exhaled air	AJKOTEKTOR PRO-100 touch-M	1 PC.
Mouthpiece	–	5 pieces. **
Funnel mouthpiece	–	1 PC. **
AA batteries 1.5 V	–	4 things.
USB cable	–	1 PC.
Case	–	1 PC.
Case	–	1 PC.
Thermal printer*	–	1 PC.
Battery pack*	–	1 PC. **
Thermal paper*	–	2 pcs. **
Charger AC adapter*	–	1 PC.
On-board charger power adapter*	–	1 PC.
Manual	RE 1.2	1 copy
Passport	–	1 copy
Notes: 1 The number of components in the table is indicated for the standard delivery set.		

End of table 6

2 Components marked with "*" are supplied when ordering analyzers complete with a thermal printer.

3 The number of components marked with "**" may differ from those indicated in the table as agreed with the Customer.

4 The delivery set includes mouthpieces supplied by the manufacturer and/or mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes."

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1.3.4 Additionally, the user is given access to download the program "Statistics PRO -100 touch - M" for internal use. The program is external software that is not metrologically significant. This program cannot lead to distortions in the analyzer measurement results displayed on the screen or printed on paper, since it is designed to collect, save and print data from the analyzer's memory on a PC and is not used by the analyzer when performing express measurements of the mass concentration of ethanol vapor in a selected sample exhaled air.

The software "Statistics PRO -100 touch - M" complies with IEC 62304 (Class A).

The program "Statistics PRO -100 touch - M", as well as instructions for working with this program are posted on the website www.alcotektor.ru in the "Software" section.

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1.3.6 Consumable materials are shown in Table 7.

Table 7 - Consumables

Name	Application	Replacement
Mouthpiece	Installed in the analyzer in measurement mode (for each subject)	Replaced after completion of the measurement procedure
Funnel mouth-piece	Installed in the analyzer in screening mode	Replaced for cleaning and disinfection and if damaged
Batteries	Analyzer internal power supply	Replaced when the resource is exhausted
Battery pack	Thermal Printer Internal Power Supply	Replaced when the resource is exhausted
Thermal paper	For thermal printer	Replaced when used

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1.4.2 The analyzer's built-in microcontroller controls the entire measurement process. The analyzer has built-in TouchM software.

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Firmware identification data is shown in Table 8.

Table 8 - Firmware Identification Data

Identification data (signs)	Meaning
Software identification name	TouchM.RU
Software version number (identification number), not lower	RU V1. 12b
Digital software identifier	F8C7DF98B38F8BBBC3 821FC6C4DCE0B9
Algorithm for obtaining a digital identifier	MD5
Note – The digital software identifier value shown in the table applies only to the firmware file of the specified version.	

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1.4.6 The analyzer's non-volatile memory allows you to store service settings data, as well as the results of the last 10,000 measurements.

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1.4.7 The analyzer blocks operation if the temperature of the measuring board on which the electrochemical sensor is installed is less than 0 °C or more plus 40 °C.

1.4.8 The analyzer has the ability to block the operation of the analyzer after the verification period has expired; in this case, when the next verification date arrives, the measurement and screening modes are blocked; For further operation, the analyzer must be submitted for periodic verification.

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1.4.9 Before sale, the analyzer is equipped with one of two variants of factory settings with a specific set of fields for data entry (in accordance with 1.4.10 or 1.4.11 of this OM, depending on the area of application of the analyzer), or an individual set of fields (in accordance with 1.4.12 of this OM).

1.4.10 The analyzers that will be used to test for alcohol intoxication by the traffic police have 6 fields for data entry with the following names (factory settings - option 1):

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1.4.11 The analyzers that will be used for medical examination for alcohol intoxication have 6 data entry fields with the following names (factory settings - option 2):

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1.4.12 For the use of analyzers in other areas (for example, for pre-trip inspections), as well as when relevant legislative regulations appear, the number and name of fields can be changed by the supplier before sale or in service centers during the operation of the analyzer.

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1.4.13 The data entered before the measurement is saved in the analyzer's memory and printed in the measurement report along with the measurement result and other data in accordance with 1.4.6 of this manual .

1.4.14 The general view of the analyzer is shown in Figure 1.

Рисунок 1 – General view of the analyzer

1.4.15 The **service button** is intended to control the operation of the analyzer in accordance with the instructions of this manual (for example, moving the cursor through the analyzer menu, entering the touch screen calibration mode, etc.).

1.4.16 The **on/off button** is designed to turn the analyzer on and off, and is also used to control the operation of the analyzer in accordance with the instructions of this manual.

1.4.17 The **touch screen** is designed to provide information to the user, to control the operation of the analyzer, as well as to enter data in measurement mode with data entry.

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1.4.18 **Connector _mini USB**, designed for data exchange between the analyzer and an external device (including a PC) using a USB cable or interface cable, as well as for connecting the analyzer power adapter (network or on-board) or the analyzer power cable.

1.4.19 Batteries are installed in the **power compartment**, which is located under a cover on the rear panel of the analyzer.

1.4.20 **mouthpiece pusher** is designed to remove the installed mouthpiece.

1.4.21 The **stylus** is inserted into the top of the analyzer body and is designed to touch the touch screen when entering data and controlling the analyzer.

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1.4.22 Between the guides for installing the mouthpiece (formed by the front and rear parts of the body), the analyzer has **two input ports of the intake system** (hereinafter referred to as fittings); through one air enters the pressure sensor, through the other air enters the electrochemical sensor (Figure 1, top view).

1.4.23 The **mouthpiece** is designed to form a flow of exhaled air entering the analyzer intake system.

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1.4.24  The **funnel mouthpiece** (Figure 3) is designed to form a flow of exhaled air for non-contact sampling.

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Рисунок 3 – Funnel mouthpiece

1.4.25 The **thermal printer** is designed for printing measurement reports on paper (thermal paper). The thermal printer has a power compartment for installing a **battery pack** and a compartment for installing a roll of **thermal paper**.

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Рисунок 4 – Appearance of the thermal printer

1.4.26 During operation, it is allowed to use the network and on-board power adapters included in the delivery kit of the analyzer, or adapters with the output characteristics specified in Table 3 and Table 4, which comply with the mandatory requirements provided for by the legislation of the Russian Federation.

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1.4.27 **USB cable** is designed to connect the analyzer to a PC when working with the “Statistics PRO -100 touch - M” program, and the USB cable is a cable no more than 3 m long with a USB plug and a mini USB plug (for connection to the analyzer).

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1.4.28 **The analyzer power cable** is designed to provide power to the analyzer from a PC and is a cable no more than 3 m long with a USB plug and a mini USB plug (for connecting to the analyzer).

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1.4.29 **The interface cable** is designed to control the analyzer and (or) ensure the transfer of measurement results to an external device.

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Depending on the model, the interface cable provides the functions listed in Table 9.

Table 9 – Functional features of interface cable models

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1.5 Marking and sealing

1.5.1 The analyzer is marked as follows:

- button designation;
- abbreviated designation of the full name of the medical product in the form of an inscription on the front panel: **ALKOTEKTOR® PRO - 100 touch - M**;
- name of the analyzer (on the rear panel): **Breath Alcohol Analyzer AJKOTEKTOR PRO-100 touch-M**;
- number of the registration certificate for the medical device;
- **SN** serial number of the analyzer;
- month and year of manufacture;
- name and contact details of the authorized representative of the manufacturer (manufacturer);
- name of the manufacturer, country of origin and contact details of the manufacturer;
- type approval sign : ;
- symbol of product type B : ;
- symbol of classification according to the degree of protection against the penetration of water and solid particles: IP20;
- sign of non-ionizing radiation : ;

– symbol “Refer to instructions for use” ;

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1.5.4 The analyzer is sealed on the left upper mounting screw in the analyzer power compartment. An authorized representative of the manufacturer carries out sealing with a label with the inscription “Breach of the seal voids the warranty”, which self-destructs when an attempt is made to open the analyzer.

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2.1.6 Do not allow condensation to accumulate in the funnel mouthpiece when using the analyzer in screening mode to prevent moisture from entering the analyzer intake system.

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2.1.10 To avoid battery failure and damage to the analyzer, **do not connect any power adapter, power cable, interface cable, or USB cable to the analyzer while batteries are installed in the analyzer's power compartment**. To operate from any of the power adapters (or cables), you must remove the batteries from the analyzer's power compartment.

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– at least two minutes must pass after smoking;

– it is recommended to rinse your mouth with water if you had a meal immediately before the measurement.

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- no interference with data transmission or control via Bluetooth ; Wireless networks operating at a frequency of 2.4 GHz may interfere with data transmission or control via Bluetooth ; radiotelephone base stations, microwave ovens and other electrical devices operating at a frequency of 2.4 GHz; metal partitions (objects) located between the thermal printer (or external device) and the analyzer;

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2.1.17 To ensure safety, if it is necessary to work with the "Statistics PRO -100 touch -M" program, the analyzer should be connected to a PC that meets the requirements of GOST IEC 60950-1.

2.2 Turning on the analyzer

2.2.1 The analyzer is turned on by pressing the button on/off, holding it until the sound signal (about 3 seconds). In this case, a splash screen appears on the screen with the abbreviated designation of the analyzer and the software version number, autotesting is carried out , after which the main menu of the analyzer appears on the screen (Figure 5).

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2.2.3 If the temperature of the measuring board is more than plus 40 ° C (less than 0 ° C), then when the analyzer is turned on, the message " **Temperature very high (low)** " appears on the screen and the analyzer turns off.

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2.3 Analyzer menu

2.3.1 The main menu of the analyzer consists of 4 icons and indicators that inform about the supply voltage level, current time and date, connection with a thermal printer (or external device) via Bluetooth, the main menu also displays the date of the next verification and information about memory fullness (Figure 5).

2.3.1.1 In the main menu of the analyzer, in the upper right corner of the screen, the current day of the week, time and date set in the analyzer are displayed in the following format:

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2.3.1.2 In the main menu of the analyzer, at the top of the screen, there is a **power indicator** , which informs about the supply voltage level. As the voltage level (discharge) of the batteries decreases, the shaded part of the indicator will decrease.


Figure 5 – Main menu of the analyzer

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2.3.3 **The screening mode** is intended for preliminary assessment of the presence of ethanol in exhaled air. In this mode, the result is given in the form of a message about the presence or absence of alcohol in the selected sample.

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-  - transition to the next menu window (when reaching the last menu window - exit to the main menu);
-  - return to the previous window;
-  - go to the main menu;
-  - printing the measurement protocol (transferring the measurement protocol to an external device);
-  - transition to the next measurement (without exiting to the main menu).

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2.5 Settings menu

2.5.1 During operation, the user has the opportunity to adjust the current date and time, select a measurement mode (with or without data entry), and also set at his own discretion: the automatic shutdown time of the analyzer, the number of copies of the printed measurement protocol, enable/disable the automatic shutdown function of the analyzer and activate the function of automatic printing of measurement protocols.

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2.5.2 To enter the settings menu in the main menu of the analyzer

(Figure 5), select the icon item . The "SETTINGS" menu will appear on the touch screen, consisting of three items intended to go to the corresponding menus (Figure 6).


Figure 6 – Settings menu

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2.5.5.4 To select the measurement mode, enter the “User Settings” menu (according to 2.5.4 of this manual) and turn the keyboard on or off by clicking on the square next to the inscription “Keyboard”, which indicates its current state:

- ☒ – Measurement mode with data input is selected
(virtual keyboard enabled)
- ☐ – Measurement mode without data input is selected
(virtual keyboard disabled)

2.5.6 Auto Shutdown Settings

To configure automatic shutdown, enter the “User Settings” menu (according to 2.5.4 of this manual) and select the “OFF” item. (Figure 7).

8). window appears on the touch screen . AUTO SHUT OFF "(Figure

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The analyzer will automatically turn off a set number of minutes after the last user action. 10 seconds before auto shutdown, “Shutdown” will appear on the analyzer screen, then a ten-second countdown will occur, accompanied by sound signals, after which the analyzer will turn off. To prevent auto shutdown at the time of counting, touch the touch screen (or press the service button), the analyzer will go to the main menu (without saving the data entered before measurement).

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2.5.6.2 If necessary, the user has the ability to disable the automatic shutdown function of the analyzer. To do this, click on the square next to the words “Auto power off”, which indicates the current state:

- ☒ Auto power off
Auto shut off function enabled
(the analyzer will turn off after the set number of minutes)
- ☐ Auto power off
Auto shut off feature disabled
(the analyzer will not turn off after the set number of minutes)

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
It is possible to set a different number of copies of the printed measurement protocol for two cases: with readings of 0.00 mg/l, and with readings of more than 0.00 mg/l (Figure 9).

2.5.7.3 Go to the “User Settings” menu (according to 2.5.4 of this OM), then click on the “PRINT” button (Figure 7), the “PRINT” window will appear on the screen (Figure 9).

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2.5.8 Bluetooth connection

2.5.8.1 Go to the “User Settings” menu (according to 2.5.4 of this manual), then click on the “Bluetooth” button (Figure 6). The “Bluetooth” window will appear on the screen (Figure 10).

2.5.8.2 In the “Bluetooth” window contains information about the Bluetooth connection status:

- “**Connection: Established**” – if there is a connection with a thermal printer or external device via Bluetooth (Figure 10, a);

- “**Connection: No**” – if there is no connection with a thermal printer or external device via Bluetooth (Figure 10, b).

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2.5.8.3 To establish a connection (pairing) of the analyzer with a thermal printer, turn on the thermal printer by long (for 2-3 seconds) pressing the on / off button on the thermal printer, and the on / off button indicator light will light up. In the “Bluetooth” window, click the “**Connect**” button, after a few seconds the connection status will change to “**Established**” (Figure 10, a), otherwise click the “**Connect**” button again.

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2.5.9 Adjusting the date and time

To adjust the current date and time in the analyzer, enter the settings menu (according to 2.5.2 of this manual). Click on the “Date and Time” button (Figure 6) – the “Date, Time Settings” window will appear on the screen (Figure 11).

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To change the date, click on the "+" or "-" buttons next to the parameter you want to change. The current time is set in a similar way.


Figure 11 – Setting the date and time

2.5.10 Engineer menu used for service settings, as well as when checking the analyzer. Entry to the menu is password protected.

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A description of how to configure the connection of the analyzer to external devices via Bluetooth or an interface cable is given in the document "Analyzers of ethanol vapors in exhaled air "ALKOTEKTOR PRO-100 touch - M" with accessories. Connecting the analyzer to external devices. Data transfer protocols ." The document is provided upon official request.

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- check the presence of all fastening elements and seals, as well as inscriptions and markings on the analyzer body - they must be clear and comply with the manual;
- make sure there are no mechanical damages that affect the performance and safety of the medical device, including the absence of external damage to the power adapters, thermal printer, charger, USB cable, analyzer power cable, interface cable;
- check the voltage level of the batteries (in accordance with 3.4.2.1 of this OM) and, if necessary, replace/charge them (in accordance with 3.4.2.2 of this OM); or connect the analyzer to the power supply (according to 1.2.16 of this manual);

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- check the correctness of the current date and time that are set in the analyzer and displayed in the main menu; if necessary, adjust the date and time in accordance with 2.5.9 of this OM.

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- check and, if necessary, adjust the analyzer's print settings in accordance with 2.5.7 of this manual;
- prepare the thermal printer for operation in accordance with 2.7.2 of this OM.

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– make sure that the analyzer is located at a distance of no more than 10 m from the external device, and there is no interference with the transmission of the measurement protocol (according to 2.1.16 of this manual);

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– if necessary, charge the battery pack (in accordance with 3.4.3 of this OM) or, if applicable, connect the thermal printer to the power supply in accordance with 1.2.17.1 of this OM.

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 *If the thermal printer is not used for more than 5 minutes, it automatically turns off in order to save battery power, and the on / off*

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
A data entry line and a virtual keyboard for typing characters will appear on the screen (Figure 12).

Note – If the analyzer's verification results have expired, and the verification work blocking function is activated in the analyzer (in accordance with 1.4.8 of this OM), then when selecting the measurement mode with data entry, the message "Requires Verification" is displayed on the screen; For further operation, the analyzer must be submitted for verification.

2.8.5 Enter data using the virtual keyboard.

When you click on the virtual button "АБВ" (or "ABC"), the virtual keyboard layout changes from Cyrillic to Latin with numbers and back (Figure 12).

Cyrillic

Latin alphabet, numbers

Figure 12 – Data entry windows in measurement mode with data entry (input field for the name of the subject)

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2.8.6 To move to the next field, press virtual

button , and if you leave the fields blank (for example, if there is no data to enter), then you can enter the missing data into the printed measurement report by hand.

Note – The entered data of fields 2, 4, 5 and 6 – option 1 of factory settings (according to 1.4.10 of this OM) or fields 3, 5 and 6 – option 2 of factory settings (according to 1.4.11 of this OM) are stored in the analyzer's memory, and when entering data before measurement, the previously entered field value is immediately displayed in the input line for these fields. Previously saved data can be corrected by deleting the previous value and entering a new one, which will also be stored until the next change.

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 **If you return to the previous field, the entered data will be lost, except for those data stored in memory.**

2.8.8 When the data in all fields has been entered, a window with the message "Insert a clean mouthpiece" will appear on the screen (Figure 13, left).


Figure 13 – Beginning of the measurement procedure

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⚠ *If the original packaging is damaged, do not use the mouthpiece.*

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When you press the “Start” button, an air sample is taken to automatically check for the absence of ethanol and other components that may affect the analyzer readings in the ambient air, in the analyzer intake system and in the mouthpiece (Figure 13, right):

– if the test result is positive (in the absence of ethanol and other components), the analyzer will go to the window ready for sampling (Figure 14 on the right);

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The exhalation should be calm and uniform (without forcing) and last about 4-5 seconds (until the end of the reference sound signal). If the

subject provides the required flow of exhaled air (at least 20 l/min), a reference sound signal sounds during exhalation, and the remaining exhalation time as a percentage is shown on the exhalation duration indicator (Figure 15). At the end of the reference sound signal, the analyzer will automatically take a sample of exhaled air, after which it is advisable to stop exhaling.

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
Figure 15 – Exhalation indication

2.8.14 If the exhaled air sample does not meet the required parameters for flow or volume (specified in 1.2.6 of this OM), i.e. in case of interruption of exhalation, the message “**Exhalation Interrupted**” is displayed on the screen. In this case, it is advisable to repeat the exhalation in accordance with 2.8.22 of this OM.

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2.8.15 If the examined person, for physiological reasons, is unable to exhale with the required parameters for flow or volume (specified in 1.2.6 of this OM), it is possible to use manual sampling in accordance with 2.9 of this OM.

2.8.16 After taking an air sample, the mass concentration of ethanol vapor in the selected sample is measured, and the message “**Analysis in progress, please wait...**” is displayed on the screen (Figure 16, left).

After the sample has been analyzed, the result of measuring the mass concentration of ethanol (alcohol) vapor in the selected sample of exhaled air is presented on the screen (Figure 16, right).


Figure 16 – Presentation of measurement results

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2.8.17 If you need to print a measurement protocol (or transfer it to an external device), click on the virtual button

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
2.8.18 Remove the used mouthpiece and dispose of it in accordance with the instructions of 4.3.2 of this manual.

2.8.19 To carry out the next measurement, press the virtual button and repeat the procedure according to 2.8.5-2.8.18 of this manual.

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2.8.20 If, when performing 2.8.11 **Alcohol detected** ” will appear on the screen (Figure 17, right).

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- return to the beginning of the measurement procedure by pressing the virtual button , in this case all data entered before the measurement will not be lost;

- make sure that the operating rules according to 2.1.4 of this manual are observed, since ethanol vapor could enter the mouthpiece from the ambient air;

- repeat the procedure according to 2.8.11 of this OM;

Alcohol detected ” appears repeatedly, remove the mouthpiece and repeat the procedure according to 2.8.9-2.8.11 of this manual.

 *As long as ethanol vapor is detected in the analyzer's mouthpiece or intake system, the analyzer will not be ready to sample.*

2.8.21 If, when performing 2.8.12 **Blow** ” window of readiness for sampling (Figure 18, left) will be replaced by the “ **No exhalation** ” window (Figure 18, right) after 30 seconds.


Figure 18 - Changing windows in the absence of exhalation

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2.8.22 If, when performing 2.8.13 of this OM, the exhaled air sample does not meet the required parameters for flow or volume (specified in 1.2.6 of this OM), i.e. in case of interruption of exhalation, the “ **Exhalation Interrupted** ” window will be displayed on the screen (Figure 19, right).


Figure 19 – Breathing test failure

It is advisable to stop exhaling. To resume the measurement procedure, click on the “ **Repeat** ” virtual button ; after the analyzer enters the mode of readiness for sampling

(Figure 1-4, right), you can repeat the exhalation in accordance with 2.8.13 of this manual.

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To do this, when the “**Exhalation Interrupted**” window is displayed on the screen (Figure 19, right), press the virtual button “**Exhalation Interrupted**” - a corresponding message will appear on the screen (Figure 20).


Figure 20 – Result “Exhalation Interrupted”

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To do this, during the exhalation of the subject, press the virtual button “**MANUAL SAMPLING**” (Figure 15), and the analyzer will perform a manual sampling. Perform manual sampling at the very end of the exhalation that the subject is capable of.

2.9.2 The successful use of exhaled air as a diagnostic medium for determining the alcohol content in a person’s blood is made possible by the fact that the amount of alcohol evaporating into the exhaled air depends on

its content in the venous blood. This dependence is a consequence of Henry’s law, which states that the pressure of a volatile (gaseous) component in a solution and, therefore, the vapor content of this component above the solution is directly proportional to its content in the solution.

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In manual sampling mode, the exhalation parameters specified in 1.2.6 of this OM are not provided.

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2.10 Refusal to measure

2.10.1 **Blow**” window of readiness for sampling will be replaced by the “**No exhalation**” window (Figure 21, left).

2.10.2 If the subject refuses to take measurements (the subject refuses to exhale), the analyzer provides the ability to interrupt the measurement procedure, while recording the fact of refusal to take measurements (if necessary). To do this, click on the virtual button “**Refusal**” - a message about refusal of measurements will appear on the screen (Figure 21, right).

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
Figure 21 – Result of “Refusal of measurements”

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2.11.1 Before taking measurements, make sure that the conditions listed in 2.1 of this OM are met, prepare the analyzer for operation in accordance with 2.6.2 of this OM, familiarize the subject with the procedure for taking measurements, as well as the rules for exhalation (in accordance with the instructions given in 2.8.13 of this RE).

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In the main menu, select the measurement mode by clicking on the corresponding icon . A window with the message “Insert a clean mouthpiece” will appear on the screen (Figure 22, left).


Figure 22 – Starting the measurement procedure in measurement mode without entering data

Note – If the analyzer's verification results have expired, and the verification work blocking function is activated in the analyzer (in accordance with Notes 1.4.8 of this OM), then when selecting the measurement mode without entering data, the message “Requires Verification” is displayed on the screen, for further operation the analyzer must be submitted for periodic verification.

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2.12.3 Select screening mode .

A message indicating that you are ready for the test will appear on the touch screen (Figure 23, left).


Figure 23 – Screening mode

Note – If the analyzer's verification results have expired, and the verification work blocking function is activated in the analyzer (in accordance with notes 1.4.8 of this OM), then when selecting the screening mode, the message "Requires Verification" is displayed on the screen, for further operation, the analyzer must be submitted for verification.

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2.12.6 If the subject provides the required flow of exhaled air, a reference sound signal sounds during exhalation, and the remaining exhalation time as a percentage is shown on the exhalation duration indicator (Figure 23, right). It is advisable to stop exhalation at the end of the reference sound signal. The analyzer will automatically take a sample of exhaled air from the funnel mouthpiece.

 ***Do not remove the analyzer from the subject's mouth until the sample is taken. Otherwise, the ambient air will be sampled rather than the exhaled air.***

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2.12.8 Depending on the test result, one of two messages appears on the analyzer screen (Figure 24).

 ***The screening mode can only be used to assess the presence of ethanol in exhaled air, and the analyzer error may exceed the limits of the permissible analyzer error specified in 1.2.1 of this OM.***


Figure 24 – Test results in screening mode

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Record” window will appear on the touch screen with the protocol of the last measurement or the last viewed protocol (Figure 25).


Figure 25 – Viewing test memory

2.13.3 If measurements were made in measurement mode with data entry, the measurement results are stored in memory together with automatically saved data (in accordance with 1.4.6 of this OM) and with data entered into the fields intended for filling in (Figure 25, left).

2.13.4 If measurements were made in the measurement mode without entering data, only the measurement results are stored in memory together with automatically saved data (according to 1.4.6 of this manual) (Figure 25, right).

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Note – Viewing records in the analyzer's memory (through one or through 10 tests) is carried out by the numbers of measurement protocols without taking into account the date and time of their implementation. It should be taken into account that after the analyzer has carried out more than 10,000 measurements (tests), subsequent protocols are recorded in the analyzer's memory instead of those saved, starting with number 00001. For example, it may turn out that after the measurement protocol with number 00001 with the date of December 15, 2021, the next In the analyzer memory there will be a measurement protocol with number 00002 with a date of October 10, 2020 (protocol data with number 00002 will be replaced with new ones only after another measurement).

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– while in the “Record” window (Figure 25), press and hold the service button and briefly press the on/off button, release both buttons; a window for switching to the required protocol will appear on the screen (Figure 26);

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– to quickly go to the measurement protocol with a specific number, enter the required test number using the “+” or “-” buttons (the most significant digit of the number is indicated on the left) and click on the virtual “Go” button. If the test with the dialed number does not exist, the transition is made to the last measurement protocol in memory.


Figure 26 – Quick transition to the required protocol

2.13.8 In the test memory viewing menu, depending on the connection status of the analyzer with a thermal printer (or an external device via Bluetooth), the title of the “**Recording**” window displays the words “**Connected**” or “**Disconnected**”:

Connected” status – the analyzer is ready to print or transfer the measurement protocol to an external device via Bluetooth;

Disabled” status – the analyzer is not ready to print or transfer the measurement protocol to an external device via Bluetooth.

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2.14.1 In the analyzer measurement protocols, information is printed according to Table 10 (an example is given with a set of fields for data entry in accordance with 1.4.10 of this OM).

Table 10 – Contents of the measurement protocol

Inscription in the protocol	Contents of the protocol
АЛКОТЕКТОР PRO-100 touch-M	Analyzer type designation
Device Number:	Analyzer serial number
Adjustment date DD/MM/YYYY	Date of last adjustment of analyzer readings (day/month/year)
Verification date: DD/MM/YYYY	Date of the last analyzer verification (day/month/year)
Test NO.:	Measurement number (according to the internal numbering of the analyzer)
Date: DD/MM/YYYY	Date of measurement (day/month/year)
Time: HH:MM	Measurement time (hour:minutes)
Mode: Automatic	Air sampling mode ¹⁾
Result: X.XX mg/l	The result of measuring the mass concentration of ethanol vapor in a selected sample of exhaled air: numerical value and designation of the unit of measurement “mg/l” ²⁾
Subject's name: ⁵⁾	Last name and initials of the person being examined ³⁾
Survey Location: ⁵⁾	Measuring location data ³⁾
State Vehicle Number: ⁵⁾	State number of the vehicle ³⁾
Badge: ⁵⁾	Inspector badge number ³⁾
Inspector: ⁵⁾	Inspector's surname and initials ³⁾
Traffic police department: ⁵⁾	Traffic police department number ³⁾
Signature Survey :	Signature of the person being examined ⁴⁾
Signature:	Inspector's signature ⁴⁾

End of table 10

- ¹⁾ In the manual air sampling mode, the message "Mode: MANUAL SAMPLING" is displayed.
- ²⁾ If it is recorded that the subject refuses to take the measurement, the message "Refusal of the test" is displayed; in the event of a recorded fact of insufficient flow and exhalation volume, the message "Exhalation Interrupted" is displayed. In this case, information about the air sampling mode is not displayed in the protocol.
- ³⁾ Data is entered from the analyzer's virtual keyboard (touch screen) before measurement or entered by hand into a printed measurement protocol.
- ⁴⁾ The data is entered by hand into the printed measurement protocol.
- ⁵⁾ The set of fields for data entry may differ from those indicated in the table (the number of fields is set from 0 to 10, the names of the fields can be changed). The names of the fields for entering measurement protocol data are indicated in the analyzer passport.

2.14.2 When transmitting a measurement result to an external device via Bluetooth or an interface cable, the transmitted protocol contains information according to Table 10 (except for signature fields - these fields are not transmitted).

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2.17 Possible malfunctions of the analyzer and thermal printer

Possible malfunctions of the analyzer and thermal printer, as well as methods for eliminating them, are given in Table 11.

Table 11 - Possible faults

No.	Possible malfunction	Probable cause	Way elimination
1	The analyzer does not turn on	1. No power supply	Install/charge batteries in accordance with 3.4.2 of this OM or connect the analyzer to the power supply in accordance with 1.2.16 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
2	The real time clock has stopped (the current date and time are not set)	1. Lithium battery undervoltage 2. The analyzer is faulty	The analyzer must be sent to a service center*
3	The analyzer does not automatically take a sample	1. Insufficient flow or volume of exhaled air	Repeat exhalation in accordance with 2.8.13 or 2.12.5 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
4	There is no sound signaling of the analyzer operation	The analyzer is faulty	The analyzer must be sent to a service center *

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End of table 11

9	When turning on the analyzer, the message "Temperature very high (low)" appears and the analyzer turns off	1. Temperature of the measuring board is more than + 40 °C (less than 0 °C)	Maintain the analyzer under operating conditions (specified in 1.1.6 of this manual) and turn it on again
		2. The analyzer is faulty	The analyzer must be sent to a service center *
10	After pressing the "Start" button, the analyzer does not enter the window ready for sample collection; the message "Alcohol detected" appears on the screen.	1. Presence of ethanol in the mouthpiece or intake system of the analyzer	Act in accordance with 2.8.20 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
eleven	The analyzer does not respond correctly or does not respond to pressing virtual buttons on the touch screen	1. Touch screen calibration reset	Calibrate the touch screen in accordance with 3.4.5 of this manual
		2. The analyzer is faulty	The analyzer must be sent to a service center*

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⚠ Install batteries of the same type and manufacturer in the power compartment. When replacing batteries, replace all 4 batteries at the same time.

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b) insert a new roll of thermal paper, positioning the roll as shown in Figure 27;

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
Figure 27 – Appearance of the paper compartment cover of a thermal printer

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Note – Touching the analyzer's touch screen with the stylus should be done gently, but with sufficient force (for example, as when writing with a pen on paper).

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– checking the analyzer readings ⁸ – it is recommended before verification (or more often in accordance with the frequency set by the user);

– adjustment of analyzer readings ⁸ – in case of negative results of checking the readings or during verification as necessary.

Note - Checking and adjusting readings can be performed by authorized technical service specialists of the organization operating the analyzer, if they have the maintenance equipment specified in Appendix E (except for a multimeter).

3.5.2 A description of the analyzer maintenance operations is given in the document “Analyzers of ethanol vapor in exhaled air “ALCOTEKTOR PRO-100 touch - M” with accessories. Maintenance instructions. Engineer’s menu” , which is available to service centers upon additional request.

Note – For organizations that have their own maintenance service and maintenance facilities, upon additional request, an extract from the document “Analyzers of ethanol vapors in exhaled air “ALCOTEKTOR PRO-100 touch -M” with accessories can be provided . Maintenance instructions. Engineer’s menu”, which contains all the necessary information for checking and adjusting readings.

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3.5.4 It is recommended to mark records of performed maintenance and repairs in the analyzer passport in the maintenance accounting table, for example, with the entry “Checking readings”, “Adjusting readings”, etc.

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4 TRANSPORTATION , STORAGE AND DISPOSAL

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4.1.3 At a minimum, it is recommended that handling signs corresponding to the following meanings be applied to transport packaging: "Fragile, be careful!", "Keep away from moisture."

4.1.4 After transportation in subzero temperatures, it is recommended to keep the analyzer in operating conditions for at least an hour before turning it on. If the "Low temperature" message appears on the analyzer screen when turned on, the holding time should be increased to two hours.

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4.2.4 The battery pack must be stored in a charged state (the ideal state of charge is 40%). If a battery pack is left to be stored discharged for a period of time of 3 months or longer, an irreversible drop in its capacity may occur, and corrosion of the cells may also occur.

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5.2 Verification of the analyzer is carried out according to the document MP-242-1670-2013 "Breath Alcohol Analyzer AJKOTEKTOR PRO-100 touch-M. Verification method".

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6.1 manufacturer's warranty

ALKOTEKTOR LLC, being an authorized representative of the manufacturer, accepts claims, complaints and reviews regarding the quality of the analyzer, and also bears warranty obligations in accordance with the provisions set out in the corresponding section of the analyzer passport.

6.2 Warranty repair

Warranty repairs of analyzers are carried out at the service center of ALKOTEKTOR LLC at the address: 199178, St. Petersburg, emb. Smolenskiy River, 5-7, letter A, tel. (812) 320-22-97, email address: www.alcote.ru, email: info@alcote.ru.

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The current list of service centers with which ALCOTEKTOR LLC has entered into a cooperation agreement is posted on the website www.alcote.ru in the "Service in the Regions" section.

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APPENDIX A

Regulations

1 IEC 60601-1:2020 "Medical electrical equipment. Part 1. General requirements for basic safety and essential performance" (GOST R IEC 60601-1-2022 " Products medical electric. Part 1. General requirements security taking into account main functional characteristics ").

2 IEC 62366-1:2015 "Medical devices – Part 1: Application of usability engineering to medical devices" (GOST R IEC 62366-1-2021 " Products medical _ Part 1. Design medical products With taking into account operational suitability ").

3 IEC 62304:2006 "Medical device software – Software life cycle processes" (GOST IEC 62304-2022 " Products medical _ Software provision. Processes life cycle ").

4 ISO 14971:2019 "Medical devices – Application of risk management to medical devices" (GOST ISO 14971-2021 " Products medical _ Application of risk management to medical devices ").

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9 Federal Law of June 26, 2008 No. 102- FZ "On Ensuring the Uniformity of Measurements."

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to IEC 61000-4-8 (GOST R 50648-94)		commercial or hospital environments
NOTE: U _n – voltage level of the electrical network before the test stimulus is applied.		

Table B. 3 - Guidance and manufacturer's declaration - interference immunity

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Recommended separation distances between portable and mobile RF communications equipment and the system

NOT INTENDED TO SUPPORT VITAL FUNCTIONS

The system is intended for use in an electromagnetic environment in which radiated interference levels are controlled. The purchaser or user of the system can avoid electromagnetic interference by maintaining a minimum separation distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, taking into account the maximum output power of the communications equipment.

Maximum transmitter output power (W)	Distance depending on transmitter frequency (m)		
	150 kHz ÷ 80 MHz	80 MHz ÷ 800 MHz	800 MHz ÷ 2.5 GHz
	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0,12	0,12	0.23
0,1	0.37	0.37	0.74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

Notes

1 At frequencies of 80 and 800 MHz, a higher field strength value is used. 2 The expressions given do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people.

3 When determining the recommended spacing values d for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power P in watts specified in the transmitter manufacturer's documentation into the given expressions

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Table C.1 - Weight and dimensional characteristics

Characteristic name	Meaning
Mouthpiece	
Length, mm, no more	111
Overall dimension in the transverse direction, mm, no more	16
Weight, g, no more	5
Mouthpiece according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes"	
Length, mm, no more	111
Overall dimension in the transverse direction, mm, no more	16
Weight, g, no more	5
Funnel mouth piece	
Height, mm, no more	40
Overall dimension in the transverse direction, mm, no more	54
Weight, g, no more	5
Thermal paper	
Width, mm	56±2
External diameter, mm	23±2
Weight, g, no more	25
Battery pack	
LxWxH, mm, no more	65x44x20
Weight, g, no more	100
Analyzer power adapter	
LxWxH (body), mm, no more	95x55x90
Cable length, m, no more	3
Weight, g, no more	450
Analyzer on-board power adapter	
LxWxH (body), mm, no more	100x35x35
Cable length (for model with cable), m, no more	3
Weight, g, no more	100
Charger	
LxWxH, mm, no more	105x85x35
Weight, g, no more	100
AC adapter charger	
LxWxH (body), mm, no more	100x60x85
Cable length, m, no more	3
Weight, g, no more	300

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	Type approval mark
	Unified sign of product circulation on the market of the Eurasian Economic Union
	On/off button
	Product symbol type B
	Symbol for classification of protection against ingress of water and solid particles
	D.C.
	Alternating current
	Non-ionizing radiation sign
	Manufacturer's name and address
	Prohibition on re-use
	Not sterile
	Do not use if packaging is damaged
	Carefully! Refer to instructions for use
	Refer to instructions for use
	Date of manufacture
	Batch number
	Don't burn

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Maintenance Tools

Table E.1 – Maintenance means

Name and type of maintenance tool	Designation of a regulatory document regulating technical requirements and (or) metrological and basic technical characteristics of a maintenance device
Combination device * Testo 622	Temperature measurement range from + 10 °C to + 30 °C, limits of permissible absolute error ± 0.5 °C. Relative humidity measurement range from 30% to 80%, permissible absolute error limits $\pm 3\%$. Absolute pressure measurement range from 80 to 110 kPa, permissible absolute error limits ± 0.5 kPa.
Generator* of gas mixtures of ethanol vapor in air GUTH model 12 V 500 complete with standard samples of the composition of aqueous solutions of ethanol VRE-2 GSO 8789-2006	Limits of permissible absolute error in reproducing the mass concentration of ethanol in gas mixtures in the range (40 – 80) mg/m ³ ± 4 mg/m ³ . Limits of permissible relative error in reproducing the mass concentration of ethanol in gas mixtures in the range (over 80 – 2000) mg/m ³ $\pm 5\%$.
Standard samples of the composition of ethanol/nitrogen gas mixtures in pressure cylinders GSO 10338-2013	Relative error limits at P=0.95 $\pm 2\%$.
Rotameter* RM-1.6 GUZ	The upper limit of volume flow measurements is 1.6 m ³ /h; limits of permissible relative error $\pm 2.5\%$ of the upper measurement limit.
Digital multimeter * ARRA-72	The DC voltage measurement range is from 0 to 60 V, the DC current measurement range is from 0 to 6 A. The permissible absolute error limits for DC voltage measurements are $\pm (0.005 U_{\text{measured}} + 0.02)$ V. Limits of permissible absolute error in direct current measurements $\pm (0.01 I_{\text{measurement}} + 0.002)$ A.
Diaphragm compressor with flow regulator	Productivity not less than 10 l/min.
Distilled water	GOST R 58144-2018.

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D.2 All maintenance equipment marked with "*" must be verified; standard samples of the composition of aqueous ethanol solutions must have valid passports.

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HS number (check points)	Nominal value of mass concentration of ethanol in HS, limits of permissible deviation, mg/l	Nominal value of the mass concentration of ethanol in aqueous solutions of ethanol ¹⁾ , limits of permissible deviation, mg/cm ³	Nominal value of the mass concentration of ethanol in HS in pressure cylinders ²⁾ , limits of permissible deviation, mg/l
1	0	distilled water	air
2	0.150 ± 0.015	0.386 ± 0.039	0.150 ± 0.015
3	0.475 ± 0.048	1.22 ± 0.12	0.475 ± 0.048 ³⁾
4	0.850 ± 0.085	2.19 ± 0.22	0.850 ± 0.085

3) In the case of using standard samples of the composition of gas mixtures C_2H_5OH/N_2 in pressure cylinders, it is allowed to use GS in a pressure cylinder with a mass concentration of ethanol from 0.33 to 0.52 mg/l as GS No. 3.

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Technical Data Sheet

1 Basic information

The Breath Alcohol Analyzer АЛКОТЕКТОР PRO-100 touch-M (hereinafter referred to as the analyzer) is designed for express measurements of the mass concentration of ethanol vapors in a selected sample of exhaled air.

Contraindications for use or side effects.

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Number in the Federal Information Foundation for Ensuring the Uniformity of Measurements 56632-14.

Manufacturer: Shenzhen Well Electric Co., Ltd, China.

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Authorized representative of the manufacturer (manufacturer): Limited Liability Company ALCOTEKTOR LLC (ALCOTEKTOR LLC), info @ alcotektor . ru .

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2.1 Measuring range and permissible error limits of the analyzer at ambient temperatures from plus 15 °C to plus 25 °C inclusive are given in Table 1.

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Measuring range mass concentration ethanol, mg/l	Limits of permissible error at temperatures from +15 to +25 °C incl.	
	absolute	relative
from 0 to 0.50 incl.	±0.05 mg/l	—
St. 0.50 to 0.95	—	±10 %

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2 If the analyzer readings exceed the upper reading limit, the message ">3.00 mg/l" is displayed on the analyzer screen and in the printed measurement protocol.

2.3 Price of the lowest digit of the scale, mg/l: 0.01.

2.4 The limits of permissible error of the analyzer depending on the ambient temperature are given in Table 2.

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Ambient temperature	Measuring range of ethanol mass concentration, mg/l	Limits of permissible errors ¹⁾	
		absolute	relative
from 0 °C to +5.0 °C incl.	from 0 to 0.25 incl.	±0.05 mg/l	—
	St. 0.25 to 0.95	—	±20 %
St. +5.0 °C to +15.0 °C incl.	from 0 to 0.33 incl.	±0.05 mg/l	—
	St. 0.33 to 0.95	—	±15 %
St. +15.0 °C to +25.0 °C incl.	from 0 to 0.50 incl.	±0.05 mg/l ²⁾	—
	St. 0.50 to 0.95	—	±10% ²⁾
St. +25.0 °C to +40.0 °C incl.	from 0 to 0.50 incl.	±0.05 mg/l	—
	St. 0.50 to 0.95	—	±10 %

¹⁾ The table shows the limits of permissible error of the analyzer under the operating conditions given in 2.20.

²⁾ According to table 1.

2.5 There are no additional errors from the presence of unmeasured components.

2.6 Parameters of the analyzed gas mixture when a sample is supplied to the analyzer input (automatic sampling mode):

– flow rate of the analyzed gas mixture, l/min: at least 20.

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2.7 Preparation time for work after switching on, s: no more than 5.

2.8 Measurement time after sampling, s: no more than 15.

2.9 Preparation time for work after measurement, s: no more than 20.

2.10 Operating time interval of the analyzer without adjusting the readings, months: not less than 12.

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2.11 The analyzer has built-in software TouchM.RU with the version number specified in 6.1.

2.12 A thermal paper roll is installed in the thermal printer, having the following overall dimensions: width 56±2 mm, outer diameter 23±2 mm.

2.13 Analyzer operating modes:

- measurement mode with data input;
- measurement mode without data input;
- screening mode;
- data transfer mode to PC;
- mode for adjusting and checking readings;
- dry gas test mode;
- touch screen calibration mode.

2.14 The analyzer has two modes of air sampling - automatic and manual.

2.15 The analyzer stores 10,000 results of the latest measurements in its memory.

2.16 Overall dimensions of the analyzer (LxWxH), mm, no more: 150x70x35.

2.17 Weight of the analyzer (with batteries, without thermal printer), g: no more than 280.

2.18 The analyzer can be powered electrically from:

- 4 AA batteries with a nominal voltage of 1.5 V (hereinafter referred to as batteries);
- 4 AA batteries with a voltage of 1.2 V and a capacity of at least 2200 mAh (hereinafter referred to as batteries);
- analyzer power adapter with an output voltage of ≈ 5.5 V and a current of at least 500 mA, connected to an AC mains;
- from the analyzer's on-board power adapter with an output voltage of ≈ 5 V and a current of at least 500 mA, connected to the vehicle's on-board network;
- from a USB port with a voltage of ≈ 5 V and a current of at least 500 mA via a USB cable, interface cable or analyzer power cable.

2.19 Number of measurements on the analyzer without replacing/charging batteries:

- without replacing batteries - at least 3000;
- without battery charge - at least 1500 (when using batteries with a capacity of at least 2200 mAh).

2.20 Terms of Use:

- ambient air temperature, °C : from 0 to +40;
- relative humidity ambient air, %: from 10 to 100 (without moisture condensation);
- atmospheric pressure, kPa: from 84.0 to 106.7.

2.21 Service life of the electrochemical sensor installed in the analyzer, years: at least 2.

2.22 Average service life of the analyzer, years: 5.

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3.2 The analyzer has the function of blocking the analyzer operation upon expiration of the verification period:

☐ No ☐ Yes

3.3 Upon request, the analyzer has the following functional features.

3.4.1 Provides the ability to work with the interface model: _____.

It is not allowed to connect an interface cable to the analyzer whose model does not correspond to the one specified.

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- ☐ No Yes: ☐ to a thermal printer via Bluetooth
☐ to an external device via Bluetooth
☐ to an external device via an interface cable

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- ☐ No ☐ Yes

3.4.5 Provides data transfer to the PRO-100 touch-M Statistics program via a USB cable:

- ☐ No ☐ Yes

4 Completeness

4.1 The analyzer delivery set is listed in Table 3.

Table 3 - Analyzer contents

Name	Designation	Quantity
Analyzer	АЛКОТЕКТОР PRO-100 touch-M	1 PC.
Mouthpiece	-	5 pieces. **
Funnel mouthpiece	-	1 PC. **
AA batteries 1.5 V	-	4 things.
USB cable	-	1 PC.
Case	-	1 PC.
Case	-	1 PC.
Thermal printer*	-	1 PC.
Battery pack*	-	1 PC. **
Thermal paper*	-	2 pcs. **
Charger AC adapter*	-	1 PC.
On-board charger power adapter*	-	1 PC.
Manual	RE 1.2	1 copy
Passport	-	1 copy

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End of table 3

3 The number of components marked with "***" may differ from those indicated in the table as agreed with the Customer.

4 The delivery set includes mouthpieces supplied by the manufacturer and/or mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes".

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6 Components that are accessories: USB cable, funnel mouthpiece, mouthpieces according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alcotektor mouthpiece with two holes", thermal printer, battery pack, thermal paper, charger, mains and on-board power adapters for the analyzer, mains and on-board power adapters for the charger, power cable for the analyzer, interface cable.

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5.3 Records of maintenance and repairs performed, as well as records of changes to the individual settings of the analyzer installed at the time of sale according to 3.1 and 3.2, and records of changes in the functional features of the analyzer according to 3.3, are recommended to be noted in Table 4.

Note - Entries in the passport about work performed are made only by service center specialists and verifiers; If necessary, the user is recommended to keep a separate maintenance log for the analyzer.

6 Certificate of acceptance

6.1 Incoming control

Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M serial number _____ has been accepted in accordance with the current technical documentation and is recognized as suitable for use.

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6.2 Primary verification

Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M serial number _____ verified in accordance with the document MP-242-1670-2013 "Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M. Verification methodology", approved by the State Center for Testing and Testing of the Federal State Unitary Enterprise "VNIIM named after. D.I. Mendeleev" November 20, 2013, at _____

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Verifier: _____
signature initials, surname

Verification mark _____

The verification results are confirmed by information on the verification results entered into the Federal Information Fund for Ensuring the Uniformity of Measurements (FSIS "ARSHIN").

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Warranty obligations

ALKOTEKTOR LLC bears warranty obligations in relation to the
 serial number _____ for 12 (twelve) months from the date of sale specified in 7.7. When the analyzer is sold to the Buyer by the Seller, ALKOTEKTOR LLC is responsible for the quality and good condition of the analyzer for 12 (twelve) months from the date specified in 7.7, but not more than 15 (fifteen) months from the date specified in 7.7. In this case, the warranty obligations of ALKOTEKTOR LLC are valid only if there is the seal of ALKOTEKTOR LLC in 7.7 and the Seller's seal in 7.8. The implementation of warranty obligations is possible only upon presentation of this passport with the signature of the Buyer (operating organization) in 7.9.

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7.2 If defects are discovered in the product during the warranty period, you must deliver the product to the service center of ALKOTEKTOR LLC. If the claim is found to be justified, the Buyer (operating organization) or the Seller has the right to have the defects eliminated free of charge or replaced with a similar analyzer (at the discretion of ALKOTEKTOR LLC).

7.3 Acceptance of the product under warranty is carried out at the service center of ALKOTEKTOR LLC, tel. : (812) 320-22-97, at the address: 199178, St. Petersburg, emb. R. Smolenki, 5 —7.

7.4 Procedure for accepting the product under warranty.

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7.8 Date of sale of the analyzer by the Seller:

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
2023

CONTENT

1 DESCRIPTION AND OPERATION	6
1.1 General provisions	6
1.1.1 Purpose	6
1.1.2 Scope	6
1.1.3 Designated users	6
1.1.4 Contraindications for use	6
1.1.5 Residual risk assessment	6
1.1.6 Operating conditions	6
1.1.15 Type of contact with the human body	7
1.2 Metrological and technical characteristics	7
1.2.15 Electrical power supply for analyzer	10
1.2.16 Electrical power supply for thermal printer	10
1.3 Composition of the analyzer	12
1.4 Design and operation of the analyzer	14
1.5 Marking and sealing	22
1.6 Packaging	23
2 INTENDED USE	23
2.1 Operating limitations	23
2.2 Turning on the analyzer	25
2.3 Analyzer menu	26
2.4 Functions of virtual touch screen buttons	28
2.5 Settings menu	29
2.5.5 Setting the measurement mode	30
2.5.6 Auto power off settings	31
2.5.7 Print Settings	32
2.5.8 Setting up Bluetooth connection	33
2.5.9 Positioning system menu	34
2.5.10 Adjusting the date and time	36
2.6 Preparing the analyzer for operation	37
2.6.1 Putting the analyzer into operation	37
2.6.2 Preparing the analyzer before starting	38
2.7 Preparation and operation of the thermal printer	39
2.8 Operating procedure in measurement mode with data entry	41
2.9 Manual sampling mode	49
2.10 Refusal to measure	50
2.11 Operating procedure in measurement mode without data entry	51

2.12 Operating procedure in screening mode	52
2.14 Contents of the printed measurement report	58
2.15 Turning off the analyzer	59
2.16 Turning off the thermal printer	60
2.17 Possible malfunctions of the analyzer and thermal printer	60
3 MAINTENANCE	63
3.4 On-site maintenance	64
3.4.2 Charge/replace analyzer batteries	64
3.4.3 Charging the battery pack (for thermal printer)	65
3.4.4 Refilling the thermal printer with thermal paper	65
3.4.5 Touch screen calibration	66
3.4.6 Cleaning and disinfection	67
3.5 Periodic maintenance	67
4 TRANSPORTATION, STORAGE AND DISPOSAL	69
4.1 Transportation	69
4.2 Storage	69
4.3 Disposal	70
5 CHECKING THE ANALYZER	71
6 MANUFACTURER'S WARRANTIES AND REPAIRS	72
6.1 Manufacturer's warranty	72
6.2 Warranty repairs	72
6.3 Post-warranty maintenance and repair	72
APPENDIX A Regulatory documents	73
APPENDIX B Information on electromagnetic compatibility	74
APPENDIX C Weight and dimensions	78
APPENDIX D Signs and symbols used in marking	80
APPENDIX E Maintenance Tools	82
CHANGE REGISTRATION SHEET	84

This Operating Manual (hereinafter referred to as the OM) is intended to study the design, principle of operation, metrological and technical characteristics of the Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M plus (hereinafter referred to as the analyzer) and contains information necessary for its proper operation and maintenance.

This instruction manual includes a description of all functionality of the analyzer. Certain provisions of this operating manual (their applicability) should be considered taking into account the delivery package of the analyzer and its individual features (which are installed upon request and indicated in the analyzer passport).

The portable automatic analyzer is simple and easy to use, based on modern advances in microelectronics.

The use of the analyzer ensures the reliability, dating and documentation of measurement results.

Persons who have read this manual and have undergone safety training in the prescribed manner are allowed to work with the analyzer.

Manufacturer (hereinafter referred to as the manufacturer):
Shenzhen Well Electric Co., Ltd, China.

Address : 1-3F, Building 101, No. 407, HeDongCun, HengKeng, GuanCheng Community, GuanHu Street, Longhua District, Shenzhen City, Guangdong, China, Tel: 86-755-83160728, Fax: 86-755-83160467, wellzp@well-co.com.

Authorized representative of the manufacturer (manufacturer):
Limited Liability Company "AJKOTEKTOP" (LLC "ALKOTEKTOR").
Location address: 191036, St. Petersburg, st. 1st Sovetskaya, 10, letter A, room . 2-N, tel. (812) 320-22-97, email . address: www.alcotector.ru, email mail: info@alcotector.ru.

AJKOTEKTOP® is a registered trademark of ALKOTEKTOR LLC.

Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M plus is registered by the Federal Service for Surveillance in the Sphere of Healthcare and is approved for circulation in the Russian Federation, registration certificate No. _____ dated _____.

Breath Alcohol Analyzers AJIKOTEKTOP PRO-100 touch-M plus are registered with the Federal Information Foundation for Ensuring the Uniformity of Measurements under No. _____.

Permitting documents are posted on the website www.alcotector.ru¹, archived versions of permitting documents are available upon official request.

The list of regulatory documents referenced in this OM is given in Appendix A.

This OM uses the following warning symbols:

 - indicates the need for the user to pay attention to the warnings and precautions given in this paragraph of this OM;

 - indicates the need for the user to pay attention to the rules for using mouthpieces given in this paragraph of text.

1 DESCRIPTION AND OPERATION

1.1 General provisions

1.1.1 Purpose

The analyzer is designed for rapid measurements of the mass concentration of ethanol vapor in a selected sample of exhaled air.

1.1.2 Application area

The analyzer is recommended for use in the areas of government regulation to ensure the uniformity of measurements:

- when carrying out activities in the field of healthcare ;
- while ensuring road safety .

1.1.3 Intended Users

Potential consumers (users) using the analyzer for its intended purpose may be:

- medical workers during medical examinations (pre-shift , pre-trip , post-shift , post-trip) and medical examinations for intoxication (alcohol, drugs or other toxic) ;
- officials who are granted the right of state control over the safety of traffic and operation of a vehicle of the corresponding type, as well as officials of the military automobile inspection;
- officials who are granted the right to perform work to control the internal labor regulations in the organization;
- individuals.

1.1.4 Contraindications for use

contraindications for use or side effects .

1.1.5 Residual risk assessment

The analyzer design process is carried out taking into account serviceability in accordance with IEC 62366-1 and the results of software requirements analysis in accordance with IEC 62304 . The hazards associated with the development, production and use of the analyzer have been assessed in accordance with ISO 14971. Necessary risk control measures have been identified and implemented to ensure that no new hazards are introduced. The residual total risk is acceptable.

1.1.6 Conditions of use

The analyzer is used under the following conditions:

- ambient temperature, °C: from minus 5 to plus 50;
- relative humidity of the ambient air (without moisture condensation), %: from 10 to 100;
- atmospheric pressure, kPa: from 84.0 to 106.7.

¹To find the required information on the site, it is advisable to use the site map.

 **Note** – To ensure safety when operating the analyzer with a power adapter connected, the requirements for the operating ambient temperature range of the power adapter must also be taken into account.

1.1.7 Abbreviated designation of the analyzer: ALKOTEKTOR PRO-100 touch - M plus.

Note – On the front panel of the analyzer there is an abbreviated designation of the medical device in the form of an inscription "АЛКОТЕКТОР® PRO-100 touch - M".

1.1.8 The analyzer is resistant to the effects of ambient temperature and humidity, the nominal values of which are specified in 1.1.6 of this manual.

1.1.9 In terms of resistance to mechanical stress, the analyzer is classified as a wearable product that is not intended to be used while being carried.

1.1.10 In terms of safety, the analyzer complies with the requirements of IEC 60601-1 and is designed as a product with an internal power supply, type B working part.

1.1.11 Information on the electromagnetic compatibility of the analyzer is given in Appendix B.

1.1.12 The analyzer housing is classified as IP 20 in terms of protection against the penetration of water and solid particles.

1.1.13 The analyzer is designed for continuous operation.

1.1.14 The noise level generated by the analyzer does not exceed 35 dBA.

1.1.15 Type of contact with the human body

1.1.15.1 The analyzer is not sterile, does not require sterilization before or after use, and does not have sterile components.

1.1.15.2 The analyzer does not contain medicines, as well as materials of animal and (or) human origin.

1.1.15.3 Mouthpieces included in the medical device and mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (version "Alcotektor mouthpiece with two holes"), which are accessories of the medical device, are intended for one-time use.

1.2 Metrological and technical characteristics

1.2.1 The metrological characteristics of the analyzer are given in Table 1.

Table 1 - Metrological characteristics

Measuring range of ethanol mass concentration, mg/l	Limits of permissible error	
	absolute	relative
from 0 to 0.200 incl.	±0.020 mg/l	—
St. 0.200 to 1.500	—	±10%

Note – The analyzer has a programmatically set minimum range of readings that are displayed on the analyzer screen and on paper in the form of zero readings: from 0.000 to 0.020 mg/l.

1.2.2 Indication range, mg/l: from 0.000 to 3.000.

Notes:

1 When the analyzer readings exceed the upper measurement limit, the mass concentration of ethanol in the selected exhaled air sample exceeds 1.350 mg/l.

2 When the analyzer readings exceed the upper reading limit, the message ">3,000 mg/l" is displayed on the analyzer screen and in the printed measurement protocol (or transmitted to an external device).

1.2.3 Price of the lowest digit of the scale, mg/l: 0.001.

1.2.4 additional errors from the presence of unmeasured components.

1.2.5 Parameters of the analyzed gas mixture when a sample is supplied to the analyzer input (automatic sampling mode):

– flow rate of the analyzed gas mixture, l/min: not less than 20;

– sample volume of the analyzed gas mixture, l: not less than 1.2.

1.2.6 Preparation time for work after switching on, s:

– at ambient temperature St. +15 °C to +50 °C: no more than 5;

– at ambient temperature from -5 °C to +15 °C inclusive: no more than 10.

1.2.7 Measurement time after sampling a gas mixture with a mass concentration of ethanol of 0.000 mg/l (0.480 mg/l), s:

– at ambient temperature St. +15 °C to +50 °C: no more than 5 (7);

– at ambient temperature from -5 °C to +15 °C inclusive: no more than 5 (15).

1.2.8 Time to prepare for work after analyzing a gas mixture with a mass concentration of ethanol of 0.000 mg/l (0.480 mg/l), s:

– at ambient temperature St. +15 °C

to +50 °C: no more than 5 (10);
 – at ambient temperature from -5 °C
 to +15 °C inclusive: no more than 5 (30).

1.2.9 The analyzer provides² the ability to print measurement reports on a thermal printer in a specified number of copies or transfer the measurement report to external devices at a distance of up to 10 m.

Bluetooth transmitter characteristics:

- frequency, MHz: from 2400 to 2483.5;
- modulation method: GFSK, 0.5 Gaussian;
- radiation power, dBm, no more than: 4.

1.2.10 Printing of measurement reports is carried out on thermal paper installed in a thermal printer. The overall dimensions of the thermal paper roll are given in Appendix B.

1.2.11 Analyzer operating modes:

– measurement mode with data entry – mode for measuring the mass concentration of ethanol vapor in a selected sample of exhaled air with the ability to pre-fill the fields of the measurement protocol from the virtual keyboard of the touch screen;

– measurement mode without data entry – a mode for measuring the mass concentration of ethanol vapor in a selected sample of exhaled air without the possibility of pre-filling the fields of the measurement protocol from the virtual keyboard of the touch screen, while it is possible to enter data into the printed protocol by hand;

– screening mode – mode for preliminary assessment of the presence of ethanol in exhaled air;

– mode of data transfer (saved measurement results) to a personal computer (hereinafter referred to as PC);

- mode for adjusting and checking readings;
- touch screen calibration mode.

1.2.12 Sampling modes:

– automatic – the intake system of the analyzer is activated by excess pressure created by the exhalation of the subject;

– manual – the analyzer intake system is activated by the operator by pressing a button on the analyzer's touch screen or on the keypad.

1.2.13 The analyzer's memory stores the results of the last 10,000 measurements.

² Upon request, analyzers may not be equipped with a Bluetooth module.

1.2.14 Automatic shutdown time of the analyzer (set by the user), min: 1/2/5/10.

1.2.15 Electrical power supply for the analyzer

1.2.15.1 The analyzer can be powered electrically from:

– 4 AA batteries with a nominal voltage of 1.5 V (hereinafter referred to as batteries);

– 4 AA batteries with a voltage of 1.2 V and a capacity of at least 2200 mAh (hereinafter referred to as batteries);

– power adapter for the analyzer with the characteristics specified in Table 2, connected to the AC mains;

– from the analyzer's on-board power adapter with the characteristics specified in Table 2, connected to the vehicle's on-board network;

– from a USB port with a voltage of ≈ 5 V and a current of at least 500 mA via a USB cable, interface cable or analyzer power cable.

Table 2 - Analyzer Power Adapter Specifications

Adapter type	Analyzer Power Adapter Specifications	
	parameter	meaning
Analyzer power adapter	Rated input voltage	~ 220 V
	Frequency	50 Hz
	Output voltage	≈ 5.5 V
	Output current	500 mA
	Electric shock protection class	II
	Operating ambient temperature range	10°C to 30°C
Analyzer on-board power adapter	Input voltage	≈ 12 V ³
	Output voltage	≈ 5 V

1.2.15.2 Number of measurements on the analyzer without replacing/charging batteries:

– without replacing batteries – at least 3000;

– without charging batteries – at least 1500 (when using batteries with a capacity of at least 2200 mAh).

1.2.16 Electrical power supply for thermal printer

1.2.16.1 The thermal printer can be powered electrically from:

– a battery pack with a voltage of 7.4 V and a capacity of at least 1800 mA · h;

³ The input voltage is determined by the model of the on-board adapter included in the delivery package of the analyzer; it can also be ≈ 24 V, $\approx 12-24$ V (indicated on the marking of the on-board power adapter).

– a charger power adapter ⁴ with the characteristics specified in Table 3, connected to an alternating current network;

– from the on-board power adapter of charger ⁴ with the characteristics indicated in Table 3, connected to the vehicle's on-board network.

Table 3 - Characteristics of charger power adapters

Adapter type nutrition	Charger Power Adapter Specifications	
	parameter	meaning
AC adapter charger	Rated input voltage	~ 100-240 V
	Frequency	50/60 Hz
	Output voltage	==12 V
	Output current	3 A
	Electric shock protection class	II
	Operating ambient temperature range	0°C to 40°C
On-board charger power adapter	Input voltage	==12 V ⁵
	Output voltage	==12 V

1.2.16.2 The battery pack can be charged:

– in the power compartment of a thermal printer ⁶ when connecting a network or on-board power adapter to the charger;

– in a charger with the characteristics specified in Table 4, when connecting to it a network or on-board power adapter of the charger with the characteristics specified in Table 3.

Table 4 - Charger characteristics

Characteristic name	Meaning
Input voltage	==12 V
Output voltage	==8.4 V
Output current, no more	1 A

1.2.17 Overall dimensions of the analyzer (LxWxH), mm, no more: 150x70x35.

⁴ Note – The ability to power the thermal printer from the power adapters of the charger (mains or on-board) is ensured if the thermal printer has a power connector.

⁵ The input voltage is determined by the model of the on-board adapter included in the delivery package of the analyzer; it can also be =24 V, =12-24 V (indicated on the marking of the on-board power adapter).

⁶ If the thermal printer has a power connector.

1.2.18 Weight of the analyzer (with batteries , without thermal printer), g, no more than: 280.

1.2.19 Weight and dimensional characteristics of components and accessories are given in Appendix B.

1.2.20 Service life of the electrochemical sensor installed in the analyzer, years: at least 2.

1.2.21 Operating time interval of the analyzer without adjusting the readings, months: not less than 12.

The readings are adjusted during verification if necessary.

1.2.22 Average service life of the analyzer, years: 5.

1.2.23 Average time to failure, hours: 8000.

1.3 Analyzer composition

1.3.1 Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M plus:

Compound:

1. Analyzer – 1 pc.
2. Mouthpieces – no more than 5 pcs. (if necessary)
3. AA batteries 1.5 V – 4 pcs.
4. Pouch – 1 pc.
5. Case – 1 pc.

Accessories:

1. Funnel mouthpiece – no more than 5 pcs.
2. Mouthpiece according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "ALCOTEKTOR mouthpiece with two holes" – no more than 5000 pcs.
3. USB cable – 1 pc.
4. Thermal printer – 1 pc.
5. Battery pack – no more than 2 pcs.
6. Thermal paper – no more than 20 pcs.
7. Charger – 1 pc.
8. Charger power adapter – 1 pc.
9. On-board charger power adapter – 1 pc.
10. Analyzer power adapter – 1 pc.
11. Analyzer on-board power adapter – 1 pc.
12. Analyzer power cable – 1 pc.
13. Interface cable – 1 pc.

1.3.2 The analyzer delivery set is shown in Table 5.

Table 5 – Analyzer contents

Name	Designation	Quantity
Analyzer ethanol vapor in exhaled air	АЛКОТЕКТОР PRO-100 touch-M plus	1 PC.
Mouthpiece	—	5 pieces. **
Funnel mouthpiece	—	1 PC. **
AA batteries 1.5 V	—	4 things.
USB cable	—	1 PC.
Case	—	1 PC.
Case	—	1 PC.
Thermal printer*	—	1 PC.
Battery pack*	—	1 PC. **
Thermal paper*	—	2 pcs. **
Charger AC adapter*	—	1 PC.
On-board charger power adapter*	—	1 PC.
Manual	RE 1.1	1 copy
Passport	—	1 copy

Notes:

1 The number of components in the table is indicated for the standard delivery set.

2 Components marked with “**” are supplied when ordering analyzers complete with a thermal printer.

3 The number of components marked with “***” may differ from those indicated in the table as agreed with the Customer.

4 The delivery set includes mouthpieces supplied by the manufacturer and/or mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version “Alkotektor mouthpiece with two holes.”

5 Accessories supplied upon request and not included in the standard delivery package: charger, mains and on-board analyzer power adapters, AA 1.2 V batteries and battery charger, analyzer power cable, interface cable.

1.3.3 The customer (developer organization) integrating the analyzer with external devices via Bluetooth and/or via an interface cable, upon official request, is provided with information on connecting external devices to the analyzer and recognizing data transfer protocols, as well as list of control commands.

1.3.4 Additionally, the user is given access to download the program “Statistics PRO -100 touch - M plus ” for internal use. The program is external software that is not metrologically significant. This program cannot lead to distortions in the analyzer measurement results displayed on the screen or printed on paper, since it is designed to collect, save and print data from the analyzer’s memory on a PC and is not used by the analyzer when performing express measurements of the mass concentration of ethanol vapor in a selected sample exhaled air.

Software “Statistics PRO -100 touch - M plus ” complies with IEC 62304 (Class A).

Program “Statistics PRO -100 touch - M plus ”, as well as instructions for working with this program are posted on the website www.alcotector.ru in the “Software” section.

1.3.5 Upon additional order, 4 pcs. can be included in the delivery package. AA 1.2 V batteries and battery charger.

1.3.6 Consumables are shown in Table 6.

Table 6 - Consumables

Name	Application	Replacement
Mouthpiece	Installed in the analyzer in measurement mode (for each subject)	Replaced after completion of the measurement procedure
Funnel mouthpiece	Installed in the analyzer in screening mode	Replaced for cleaning and disinfection and if damaged
Batteries	Analyzer internal power supply	Replaced when the resource is exhausted
Battery pack	Thermal Printer Internal Power Supply	Replaced when the resource is exhausted
Thermal paper	For thermal printer	Replaced when used

1.4 Design and operation of the analyzer

1.4.1 The operating principle of the analyzer is based on the use of an electrochemical sensor for measuring the mass concentration of ethanol vapor in the analyzed air sample, which represents an electrochemical cell with two platinum electrodes, on the anode of which a catalyst specific to ethanol is deposited. The distinctive properties of the electrochemical sensor are specificity towards ethanol, high sensitivity, accuracy and stability.

1.4.2 The analyzer's built-in microcontroller controls the entire measurement process. The analyzer has built-in software TouchMplus .RU.

The built-in system-application software of the analyzer (hereinafter referred to as the built-in software) was developed by the manufacturer specifically to solve the problem of measuring the mass concentration of ethanol vapor in exhaled air, as well as displaying measurement results on the screen, storing measured data and transmitting measured data to external devices. The firmware version number is indicated in the analyzer passport. Firmware identification is done by displaying the version number on the screen when the analyzer is turned on.

The influence of the built-in software on the metrological characteristics of the analyzer is taken into account when standardizing them. The analyzers have firmware protection against intentional or unintentional changes. The software protection level is "medium" in accordance with the metrology recommendation R 50.2.077-2014.

The firmware complies with IEC 62304 (Class A).

Firmware identification data is shown in Table 7.

Table 7 - Firmware Identification Data

Identification data (signs)	Meaning
Software identification name	TouchMplus.RU
Software version number (identification number)	RU V1.XXp
Digital software identifier	CA743650F391C3A13FE-AADB7330CC422
Algorithm for obtaining a digital identifier	MD5
Note – The value of the digital software identifier indicated in the table applies only to the firmware file with version number RU V1.00p.	

1.4.3 The analyzer is controlled using a touch screen and two buttons located on the front panel, and (or) control commands from an external device (via Bluetooth or using an interface cable).

1.4.4 The stages of operation and air sampling of the analyzer are accompanied by sound signals.

1.4.5 The analyzer is equipped with a real-time clock chip, which is powered when there is no power supply to the analyzer from a lithium battery installed on the main board.

1.4.6 When equipped with a built-in positioning system receiver (GLONASS/GPS), the analyzer provides the ability to determine the coordinates of the measurement location and automatically enter them into the measurement protocol.

Note – The user activates the function of automatically entering the coordinates of the measurement location into the printed measurement protocol independently in accordance with 2.5.9.6 of this OM.

1.4.7 The analyzer's non-volatile memory allows you to store service settings data, as well as the results of the last 10,000 measurements.

Along with the measurement result, automatically saved data is stored in the memory (serial number of the measurement (test), date and time of its implementation, sample collection mode, date of correction of readings and calibration of the analyzer), and data entered before the measurement (name of the subject, place of examination, etc.) .d.).

Note – If by the time the measurement is taken, the receiver of the analyzer positioning system has determined the location coordinates, then these coordinates are also automatically saved along with the measurement result.

After the analyzer memory is full, subsequent measurements will be saved, replacing the measurements with the lowest serial number, and the numbering of tests will begin again from 00001.

1.4.8 The analyzer blocks operation if the temperature of the measuring board on which the electrochemical sensor is installed is less than minus 5 °C or more plus 50 °C.

1.4.9 The analyzer has the ability to block the operation of the analyzer after the verification period has expired; in this case, when the next verification date arrives, the measurement and screening modes are blocked; For further operation, the analyzer must be submitted for periodic verification.

This function is activated only at the request of the customer (or in connection with changes in legislation) in a service center or in an organization accredited to verify analyzers.

1.4.10 Before sale, the analyzer is equipped with one of two variants of factory settings with a specific set of fields for data entry (in accordance with 1.4.11 or 1.4.12 of this OM, depending on the area of application of the analyzer), or an individual set of fields (in accordance with 1.4.13 of this OM).

1.4.11 The analyzers that will be used to test for alcohol intoxication by the traffic police have 6 fields for data entry with the following names (factory settings - option 1):

- Field 1 – **Name of the Subject** ;
- Field 2 – **Survey Location** ;
- Field 3 – **State Vehicle Number** ;

Field 4 – **Badge** (inspector);

Field 5 – **Inspector** ;

Field 6 – **Traffic Police Department** .

The entered data in fields 2, 4, 5 and 6 are stored in the analyzer's memory until the user enters new data instead (in accordance with Note 2.8.6 of this OM).

1.4.12 The analyzers that will be used for medical examination for alcohol intoxication have 6 data entry fields with the following names (factory settings - option 2):

Field 1 – **Full name of the person being examined**;

Field 2 – **Year of Birth Survey** ;

Field 3 – **Sent** (by whom sent for examination);

Field 4 – **Direction number**;

Field 5 – **Survey Location**;

Field 6 – **Full name of the Health Worker**.

The entered data in fields 3, 5 and 6 are stored in the analyzer's memory until the user enters new data instead (as specified in Note 2.8.6 of this OM).

1.4.13 For the use of analyzers in other areas (for example, for pre-trip inspections), as well as when relevant legislative regulations appear, the number and

The names of the fields can be changed by the supplier before sale, or at service centers during the operation of the analyzer.

The analyzer has the ability to set from 0 to 10 fields, as well as set (check) those fields whose entered data will be stored in the analyzer's memory until the user enters new data instead.

1.4.14 The data entered before the measurement is saved in the analyzer's memory and printed in the measurement report along with the measurement result and other data in accordance with 1.4.7 of this manual .

1.4.15 Before sale, the function of entering signature fields into the measurement protocol can be activated in the analyzer; in this case, in the printed measurement protocol, after the data entry fields, signature lines for the person being examined and the person who carried out the measurement will be printed.

1.4.16 The general view of the analyzer is shown in Figure 1.

Рисунок 1 – General view of the analyzer

1.4.17 **The service button** is intended to control the operation of the analyzer in accordance with the instructions of this manual (for example, moving the cursor through the analyzer menu, entering the touch screen calibration mode, etc.).

1.4.18 **The on/off button** is designed to turn the analyzer on and off, and is also used to control the operation of the analyzer in accordance with the instructions of this manual.

1.4.19 **The touch screen** is designed to provide information to the user, to control the operation of the analyzer, as well as to enter data in measurement mode with data entry.

Data can be entered in Cyrillic, Latin letters and numbers. Selecting menu items or entering characters is done by touching the corresponding touch screen icon **with a stylus** or finger.

1.4.20 **Connector mini USB** , designed for data exchange between

the analyzer and an external device (including a PC) using a USB cable or interface cable, as well as for connecting the analyzer power adapter (network or on-board) or the analyzer power cable.

1.4.21 The batteries are installed in the **power compartment**, which is located under the cover on the rear panel of the analyzer.

1.4.22 **mouthpiece pusher** is designed to remove the installed mouthpiece.

1.4.23 The **stylus** is inserted into the top of the analyzer body and is designed to touch the touch screen when entering data and controlling the analyzer.

The stylus has a tip made of special soft plastic that does not leave scratches on the touch screen.

1.4.24 Between the guides for installing the mouthpiece (formed by the front and rear parts of the body), the analyzer has **two input ports of the intake system** (hereinafter referred to as fittings); through one air enters the pressure sensor, through the other air enters the electrochemical sensor (Figure 1, top view).

1.4.25 The **mouthpiece** is designed to form a flow of exhaled air entering the analyzer intake system.

The special shape of the mouthpiece (Figure 2) provides excess air pressure at the inlet of the analyzer intake system during exhalation and thereby provides the ability to control the flow rate and volume of exhaled air. On the underside of the mouthpiece there are two holes for the analyzer fittings.

An arrow is applied to the mouthpiece (using injection molding), indicating the direction of the flow of exhaled air.

⚠ The mouthpiece is intended for one-time use in measurement mode with and without data input.

The analyzer uses mouthpieces supplied by the manufacturer, or mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (version "Alcotektor mouthpiece with two holes").

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1.4.26 ⓘ The **funnel mouthpiece** (Figure 3) is designed to form a flow of exhaled air for non-contact sampling.

The funnel mouthpiece has two holes on the bottom side for the analyzer fittings.

The funnel mouthpiece is intended for reusable use in screening mode and is subject to disinfection in accordance with 3.4.6 of this OM.

Рисунок 3 – Funnel mouthpiece

1.4.27 The **thermal printer** is designed for printing measurement reports on paper (thermal paper). The thermal printer has a power compartment for installing a **battery pack** and a compartment for installing a roll of **thermal paper**.

The appearance of the thermal printer is shown in Figure 4.

The thermal printer has an on/off button on the front panel, and may also have a connector for connecting a network or on-board power adapter charger.

Рисунок 4 – Appearance of the thermal printer

1.4.28 During operation, it is allowed to use the network and on-board power adapters included in the delivery kit of the analyzer, or adapters with the output characteristics specified in Table 2 and Table 3, which comply with the mandatory requirements provided for by the legislation of the Russian Federation.

The length of the power adapter cable (mains or onboard) should be no more than 3 m.

1.4.29 **USB cable** is designed to connect the analyzer to a PC when working with the program "Statistics PRO -100 touch - M" plus ", and also The USB cable is a cable no longer than 3 m with a USB plug and a mini USB plug (for connecting to the analyzer).

Note - The USB cable can be used to connect the analyzer to the on-board power adapter in the form of a cigarette lighter plug (without a cable and with a USB type A connector).

1.4.30 **The analyzer power cable** is designed to provide power to the analyzer from a PC and is a cable no more than 3 m long with a USB plug and a mini USB plug (for connecting to the analyzer).

Note - The analyzer power cable can be used to connect the analyzer to the on-board power adapter in the form of a cigarette lighter plug (without a cable and with a USB type A connector).

1.4.31 **The interface cable** is designed to control the analyzer and (or) ensure the transfer of measurement results to an external device.

The interface cable has a USB plug for connecting to a PC (or an optional mini USB or USB type -C plug for connecting to an external device) and a mini USB plug for connecting to the analyzer.

Depending on the model, the interface cable provides the functions listed in Table 8.

Table 8 – Functional features of interface cable models

Functions	Interface cable model		
	IK-1	IK-2	IK-3
Electrical power supply for the analyzer	+	+	+
Transferring data (measurement results) to an external device	+	-	+
Control from an external device	-	-	+
Automatic switching on of the analyzer when power is supplied to the mini USB connector of the analyzer	-	+	+

1.5 Marking and sealing

1.5.1 The analyzer is marked as follows:

- button designation;
- abbreviated designation of the full name of the medical product in the form of an inscription on the front panel: **ALKOTEKTOR® PRO - 100 touch - M**;
- name of the analyzer (on the rear panel): Breath Alcohol Analyzer АЛКОТЕКТОР PRO-100 touch-M plus;
- number of the registration certificate for the medical device;
- **SN** serial number of the analyzer;
- month and year of manufacture;
- name and contact details of the authorized representative of the manufacturer (manufacturer);
- name of the manufacturer, country of origin and contact details of the manufacturer;
- type approval sign:
- symbol of product type B:
- symbol of classification according to the degree of protection against the penetration of water and solid particles: IP20;
- sign of non-ionizing radiation:
- symbol "Refer to instructions for use"

mini USB connector designation : ;

– safety sign “Follow the operating instructions” (at the mini USB connector) : .

1.5.2 A label containing information for identifying the analyzer among Bluetooth devices is placed in the power compartment of the analyzer, which provides transmission of the measurement protocol to an external device via Bluetooth.

1.5.3 The list of signs and symbols used in the marking of the analyzer, its components and accessories, as well as on the transport packaging, and their explanation is given in Appendix D.

1.5.4 sealed on the left upper mounting screw in the analyzer power compartment. An authorized representative of the manufacturer carries out sealing with a label with the inscription “Breach of the seal voids the warranty”, which self-destructs when an attempt is made to open the analyzer.

1.6 Package

For transportation purposes, the complete analyzer is packaged in transport containers - corrugated cardboard boxes.

Mouthpieces packed in primary packaging made of polypropylene or polyethylene film are packaged in secondary packaging (bags made of polyethylene film).

2 INTENDED USE

2.1 Operating restrictions

2.1.1 To ensure safety requirements and prevent unforeseen risks, the user must comply with and be responsible for fulfilling the requirements of 2.1.2-2.1.17.

2.1.2 Before using the analyzer, make sure that the measurement conditions meet the requirements of 1.1.6 of this manual.

2.1.3 Modification of the analyzer is not permitted.

2.1.4 It is not allowed to store and use the analyzer in premises where alcohol-containing substances are stored in open containers, or where surfaces or equipment are treated with alcohol-containing solutions.

2.1.5 It is not allowed to store and use the analyzer in rooms with increased dust and gas contamination to avoid contamination of the analyzer intake system.

2.1.6 Do not allow condensation to accumulate in the funnel mouth-piece when using the analyzer in screening mode to prevent moisture from entering the analyzer intake system.

2.1.7 It is not allowed to connect devices and cables to the analyzer that do not meet the requirements established in this manual, as well as in the analyzer's passport.

2.1.8 Using the USB cable and the power cable of the analyzer, it is not allowed to connect any other devices to the mini USB connector of the analyzer except a PC and the on-board power adapter of the analyzer (in the form of a “cigarette lighter” plug with a USB type A connector).

2.1.9 Using an interface cable, you can connect external devices to the mini USB connector of the analyzer that provide power supply from voltage ≈ 5 V and current not less than 500 mA, in accordance with the information related to electromagnetic compatibility (Appendix B).

2.1.10 To avoid battery failure and damage to the analyzer, **do not connect any power adapter, power cable, interface cable, or USB cable to the analyzer while batteries are installed in the analyzer's power compartment.** To operate from any of the power adapters (or cables), you must remove the batteries from the analyzer's power compartment.

2.1.11 It is required to follow the rules for the safe operation of devices connected to the analyzer.

The means of electrically separating the analyzer from all poles of the power supply network when powering the analyzer (thermal printer, charger) from the mains adapter is the mains plug.

2.1.12 To exclude the possible influence of ethanol on the oral mucosa on the measurement result, at least 20 minutes must pass before the measurement after consuming alcohol-containing medications and oral sprays, as well as low-alcohol foods (fermented milk products, kvass, etc.).

2.1.13 To avoid contamination of the analyzer intake system, the analyzed air sample should not contain particles of tobacco smoke, sputum (saliva) and food debris. Therefore, before measuring:

- at least two minutes must pass after smoking;
- it is recommended to rinse your mouth with water if you had a meal immediately before the measurement.

2.1.14 Do not use the mouthpiece if its original packaging is damaged.

2.1.15 During operation, it is necessary to ensure conditions to prevent the analyzer from falling out of the user's hands to avoid damage to it. Constructing the analyzer housing from impact-resistant material and without sharp protrusions reduces the risk of damage from sharp protrusions or debris if the analyzer is dropped. For additional protection, you can use the analyzer placed in a case.

2.1.16 For working with a thermal printer (or for transmitting data via Bluetooth to external devices) conditions for normal functioning of the analyzer together with it should be ensured, namely:

- no interference with data transmission or control via Bluetooth. Wireless networks operating at a frequency of 2.4 GHz may interfere with data transmission or control via Bluetooth; radiotelephone base stations, microwave ovens and other electrical devices operating at a frequency of 2.4 GHz; metal partitions (objects) located between the thermal printer (or external device) and the analyzer;

- the thermal printer (or external device) must be located no more than 10 meters from the analyzer.

2.1.17 To ensure safety, if it is necessary to work with the "Statistics PRO -100 touch -M plus" program, the analyzer should be connected to a PC that meets the requirements of GOST IEC 60950-1.

2.2 Turning on the analyzer

2.2.1 The analyzer is turned on by pressing the button on/off, holding it until the sound signal (about 3 seconds). In this case, a splash screen appears on the screen with the abbreviated designation of the analyzer and the software version number, autotesting is carried out, after which the main menu of the analyzer appears on the screen (Figure 5).

2.2.2 When the analyzer is turned on, the date of the next verification is automatically checked. If the verification results have expired, the following message appears on the screen when the analyzer is turned on, accompanied by an audible signal:

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This message is a reminder to submit the analyzer for inspection.

2.2.3 If the temperature of the measuring board is more than plus 50 ° C (less than minus 5 ° C), then when you turn on the analyzer, the message "Temperature very high (low)" appears on the screen and the analyzer turns off.

Maintain the analyzer under the measurement conditions specified in 1.1.6 of this manual and turn it on again.

2.3 Analyzer menu

2.3.1 The main menu of the analyzer consists of 4 icons and indicators informing about the supply voltage level, current time and date, establishing the coordinates of the current location by the positioning system receiver, connecting to a thermal printer (or external device) via Bluetooth, and the date of the next verification is also displayed in the main menu and information about memory capacity (Figure 5).

2.3.1.1 In the main menu of the analyzer, in the upper right corner of the screen, the **current day of the week, time and date** set in the analyzer are displayed in the following format:

Day of week Hours :Minutes:Seconds

Day month Year

The current date and time are recorded when a sample is taken, stored in the analyzer's memory along with the measurement result and printed in the measurement report.

2.3.1.2 In the main menu of the analyzer, at the top of the screen, there is a **power indicator**, which informs about the supply voltage level. As the voltage level (discharge) of the batteries decreases, the shaded part of the indicator will decrease.


Figure 5 – Main menu of the analyzer

2.3.1.3 In the main menu of the analyzer, at the bottom of the screen, the following information is displayed:

- **date of the next verification** (the verification date is set or adjusted in the analyzer's memory if the verification results are positive);
- information about **memory capacity** (the number of measurements (tests) stored in the analyzer's memory).

2.3.1.4 **Bluetooth indicator** appears in the main menu of the analyzer, in the upper left part of the screen, when the analyzer is connected (paired) to a thermal printer or external device via Bluetooth.

2.3.1.5 **The positioning system indicator** appears in the analyzer's main menu, at the top left of the screen, when the current location coordinates are established.

2.3.2 **The measurement mode** is intended for measuring the mass concentration of ethanol vapor in a selected sample of exhaled air. Depending on the settings of the analyzer (according to 2.5.5 of this manual), when you click on the measurement mode icon, you switch to the measurement mode with data input or to the measurement mode without data input.

2.3.3 **The screening mode** is intended for preliminary assessment of the presence of ethanol in exhaled air. In this mode, the result is given in

the form of a message about the presence or absence of alcohol in the selected sample.

⚠ The screening mode can only be used to assess the presence of ethanol in exhaled air, and the analyzer error may exceed the limits of the permissible analyzer error specified in 1.2.1 of this OM.

The results of tests carried out in screening mode are stored in the analyzer's memory, but they cannot be printed on a printer in the form of protocols.

Note – If the analyzer is configured to transfer the measurement protocol to an external device via Bluetooth or when using an interface cable, the results of tests carried out in screening mode are automatically transferred to the external device after being presented on the analyzer screen.

In screening mode, it is advisable to use a reusable funnel mouthpiece included in the analyzer delivery package.

2.3.4 **memory viewing menu** allows you to view the tests recorded in the memory (saved measurement results), as well as print or transfer their protocols to an external device (according to 2.13 of this manual).

2.3.5 **The settings menu** includes a user settings menu, a date and time setting menu, and an engineer menu.

2.3.6 To select a menu item, as well as activate virtual buttons on the touch screen, touch them with a stylus (finger) or move the cursor (red frame) using the service button to the desired virtual button or menu item and confirm the selection by briefly pressing the analyzer on/off button.

2.4 Functions of virtual touch screen buttons

To navigate between menu windows, as well as print (transfer to an external device) the measurement protocol, the following virtual control buttons located at the bottom of the touch screen are used:

- transition to the next menu window (when reaching the last menu window - exit to the main menu);

- return to the previous window;

- go to the main menu;

- printing the measurement protocol (transferring the measurement protocol to an external device);

- transition to the next measurement (without exiting to the main menu).

To activate the desired button, touch it with the stylus or move the cursor (red frame) using the service button to the desired virtual button on the touch screen and confirm the selection by briefly pressing the analyzer on/off button.

2.5 Settings menu

2.5.1 During operation, the user has the opportunity to adjust the current date and time, select a measurement mode (with or without data entry), and also set at his own discretion: the automatic shutdown time of the analyzer, the number of copies of the printed measurement protocol, enable/disable the automatic shutdown function of the analyzer, activate the function of automatically printing measurement protocols and automatically entering the coordinates of the measurement location into the measurement protocol.

The above settings are made in the corresponding menus, which are entered from the settings menu.

2.5.2 To enter the settings menu in the main menu of the analyzer (Figure 5), select the icon item . The "SETTINGS" menu will appear on the touch screen, consisting of three items intended to go to the corresponding menus (Figure 6).



Figure 6 – Settings menu

2.5.3 Bluetooth connection settings, as well as entering the positioning system menu are made in the "User Settings" menu.

2.5.4 To enter the "User Settings" menu, enter the settings menu according to 2.5.2 of this OM, then click on the "User Settings" virtual button. A menu consisting of five items will appear on the touch screen (Figure 7).

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Figure 7 – User settings menu

2.5.5 Setting the measurement mode

2.5.5.1 During operation, the user has the ability to set one of two measurement modes:

- measurement mode with data input;
- measurement mode without data input.

To set the measurement mode with data entry, you must enable the virtual keyboard; for the mode without data entry, disable it (according to 2.5.5.4 of this manual).

2.5.5.2 The measurement mode with data entry allows you to carry out documented measurements. This is achieved thanks to the ability to enter data before measurement (filling in the fields of the measurement protocol) using a virtual keyboard on the touch screen, and saving the entered data in the analyzer's memory.

2.5.5.3 The No Data Input measurement mode is used to take measurements without the need to enter data before taking a measurement. In this mode, filling in the fields is only possible by hand in the printed measurement protocol.

2.5.5.4 To select the measurement mode, enter the “**User Settings**” menu (according to 2.5.4 of this manual) and turn the keyboard on or off by clicking on the square next to the inscription “**Keyboard**”, which indicates its current state:

- ☒ – **Measurement mode with data input** is selected
(virtual keyboard enabled)
- ☐ – **Measurement mode without data input** is selected
(virtual keyboard disabled)

2.5.6 Auto Shutdown Settings

To configure automatic shutdown, enter the “**User Settings**” menu (according to 2.5.4 of this manual) and select the “**OFF**” item.” (Figure 7).

7). window appears on the touch screen . **AUTO SHUT OFF**”(Figure 8).


Figure 8 – Automatic shutdown settings

2.5.6.1 To set the analyzer's automatic shutdown time, select one of the settings by touching the field to the left of the selected value.

The analyzer will automatically turn off a set number of minutes after the last user action. 10 seconds before auto shutdown, “**Shutdown**” will appear on the analyzer screen, then a ten-second countdown will occur, accompanied by sound signals, after which the analyzer will turn off. To prevent auto shutdown at the time of counting, touch the touch screen (or press the service button), the analyzer will go to the main menu (without saving the data entered before measurement).

It is recommended to select the optimal automatic shutdown time to save power.

2.5.6.2 If necessary, the user has the ability to disable the automatic shutdown function of the analyzer. To do this, click on the square next to the words “**Auto power off**”, which indicates the current state:

- ☒ **Auto power off** **Auto shut off function enabled**
(the analyzer will turn off after the set number of minutes)
- ☐ **Auto power off** **Auto shut off feature disabled**
(the analyzer will not turn off after the set number of minutes)

2.5.7 Print Settings

2.5.7.1 During operation, the user, at his discretion, can activate the function of automatic printing of measurement reports in a specified number of copies.


In this case, the thermal printer will automatically print the measurement protocol in a specified number of copies immediately after the measurement result is displayed on the screen.

If the automatic printing function is not activated, then measurement protocols will be printed only by clicking on the virtual button  at the bottom of the screen.

2.5.7.2 During operation, the user, at his own discretion (or in accordance with legal requirements), has the opportunity to set the number of copies of the printed measurement report, which will be printed automatically after the measurement result is displayed on the screen. The maximum number of protocol copies is 5.

It is possible to set a different number of copies of the printed measurement protocol for two cases: with readings of 0.000 mg/l, and with readings of more than 0.000 mg/l (Figure 9).

2.5.7.3 Go to the “**User Settings**” menu (according to 2.5.4 of this OM), then click on the “**PRINT**” button (Figure 7), the “**PRINT**” window will appear on the screen (Figure 9).

Drawing 9 – Print Settings

2.5.7.4 Set the number of copies by pressing the “+” and “-” buttons until the desired number is reached.

In order for the measurement protocol to be printed in a specified number of copies automatically after the measurement result is presented on the screen, you should activate the automatic printing function by clicking on the “**AUTO PRINT**” box (Figure 9).

2.5.8 Bluetooth connection

2.5.8.1 Go to the “**User Settings**” menu (according to 2.5.4 of this manual), then click on the “**Bluetooth**” button (Figure 6). The “**Bluetooth**” window will appear on the screen (Figure 10).

2.5.8.2 In the “**Bluetooth**” window contains information about the Bluetooth connection status :

- “**Connection: Established**” – if there is a connection with a thermal printer or external device via Bluetooth (Figure 10, a) ;

- “**Connection: No**” – if there is no connection with a thermal printer or external device via Bluetooth (Figure 10, b) .


Figure 10 – Bluetooth window (possible options)

2.5.8.3 To establish a connection (pairing) of the analyzer with a thermal printer, turn on the thermal printer by long (for 2-3 seconds) pressing the on / off button on the thermal printer, and the on / off button indicator light will light up. In the “**Bluetooth**” window, click the “**Connect**” button, after a few seconds the connection status will change to “**Established**” (Figure 10, a), otherwise click the “**Connect**” button again.

2.5.8.4 When a connection (pairing) of the analyzer with a thermal printer is successfully established via Bluetooth The connection status will change to “**Established**”, and a Bluetooth indicator will appear in the main menu of the analyzer, and in the test memory viewing menu, the title of the “**Recording**” window will show “**Connected**”.

If there is no connection between the analyzer and the thermal printer after completing 2.5.8.3 of this manual, make sure that the conditions specified in 2.1.16 of this manual are met.

2.5.9 Positioning system menu

2.5.9.1 IN In the positioning system menu, the user can view the current location coordinates (longitude and latitude) and, if necessary, activate the function of automatically entering them into the measurement protocol.

2.5.9.2 The built-in positioning system receiver (GLONASS/GPS) (hereinafter referred to as the receiver) is activated when the analyzer is turned on. If the analyzer has been turned off for the last half hour, initializing the receiver may take about 5-10 minutes.


For stable and accurate operation, the analyzer receiver needs “open sky”, that is, the absence of interference above it and nearby (buildings, trees, tunnels, etc.). The operation of the analyzer's receiver is difficult on cloudy days and in densely built areas. When using the analyzer in buildings, tunnels, etc., where the satellite signal is usually absent or extremely weak, the operation of the receiver may also be difficult.

When interference occurs, the receiver may either not determine the coordinates at all, or determine them with less accuracy. The level of accuracy depends on many factors and is not guaranteed.

The time it takes for the receiver to determine the location coordinates can be 5-10 minutes or more, depending on the above conditions.

2.5.9.3 When the location coordinates have been established by the receiver, a positioning system indicator will appear in the upper left corner of the analyzer's main menu screen.

2.5.9.4 User Settings” menu (according to 2.5.4 of this OM), then click on the “Coordinates” button (Figure 7). The “Coordinates” window will appear on the screen (Figure 11).

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Figure 11 – Positioning system menu

⚠ After the measurement is completed, the receiver stops updating the location coordinates, while the positioning system menu stores the analyzer location coordinates recorded at the time the measurement result was displayed on the screen. Therefore, **to update the location coordinates, you must turn off and turn on the analyzer again**.

2.5.9.5 If by the time the measurement is taken, the analyzer receiver has determined the location coordinates, then these coordinates are saved along with the measurement result.

2.5.9.6 To activate the function of automatically entering the coordinates of the measurement location into the measurement protocol, click on the “Print” square (Figure 11). In this case, a field with the coordinates of the measurement location will be added to the printed measurement protocol.

2.5.9.7 To display the coordinates of the measurement location on the map, you can use any free Internet service, where in the line for entering the address/coordinates, enter the coordinates specified in the printed measurement protocol and press the search button.

When entering coordinates in the search bar, use Latin letters, use spaces to separate degrees and minutes, use a dot to separate the whole and fractional parts of minutes, and enter coordinates separated by commas. An example of the format for entering coordinates is shown in Figure 12.

E030 15.69541, N59 56.94355

Figure 12 – Example of entering coordinates

2.5.10 Adjusting the date and time

To adjust the current date and time in the analyzer, enter the settings menu (according to 2.5.2 of this manual). Click on the “Date and Time” button (Figure 6) – the “Date, Time Settings” window will appear on the screen (Figure 13).

You can also enter the “Date, Time Settings” window by clicking on the area in the upper right corner of the screen in the main menu of the analyzer, where the current date and time are displayed (Figure 5).

To change the date, click on the “+” or “-” buttons next to the parameter you want to change. The current time is set in a similar way.


Figure 13 – Setting the date and time

2.5.1.1 **Engineer menu** used for service settings, as well as when checking the analyzer. Entry to the menu is password protected.

The engineer's menu includes: adjusting readings, setting parameters of the sampling system, setting a set of fields for data entry, setting the verification date in the analyzer's memory and other options.

The menu description is given in the document "Analyzers of ethanol vapors in exhaled air "ALKOTEKTOR PRO-100 touch - M" with accessories. Maintenance instructions. Engineer's menu", which is available upon additional request (in accordance with 3.5 of this OM).

2.6 Preparing the analyzer for operation

2.6.1 Commissioning the analyzer

2.6.1.1 To put the analyzer into operation, you must check the completeness of the analyzer in accordance with the passport and/or delivery agreement (if any) and check the information about the results of verification of the analyzer in the Federal Information Foundation for Ensuring the Uniformity of Measurements.

2.6.1.2 When connecting the analyzer to an external device to transfer data (measurement results) and/or control the analyzer using control commands, you must configure the connection to the external device.

A description of how to configure the connection of the analyzer to external devices via Bluetooth or an interface cable is given in the document "Analyzers of ethanol vapors in exhaled air "ALKOTEKTOR PRO-100 touch - M" with accessories. Connecting the analyzer to external devices. Data transfer protocols." The document is provided upon official request.

2.6.2 Preparing the analyzer before starting work

Before starting work, carry out an external inspection of the analyzer and prepare for work:

- check the presence of all fastening elements and seals, as well as inscriptions and markings on the analyzer body - they must be clear and comply with the manual;

- make sure there are no mechanical damages that affect the performance and safety of the medical device, including the absence of external damage to the power adapters, thermal printer, charger, USB cable, analyzer power cable, interface cable;

- check the voltage level of the batteries (in accordance with 3.4.2.1 of this OM) and, if necessary, replace/charge them (in accordance with 3.4.2.2 of this OM); or connect the analyzer to the power supply (according to 1.2.15 of this manual);

⚠ When powering the analyzer from an external source (using a power adapter, interface cable, power cable, or USB cable), be sure to remove the batteries from the analyzer's power compartment.

- check the correctness of the current date and time that are set in the analyzer and displayed in the main menu; if necessary, adjust the date and time in accordance with 2.5.10 of this OM.

⚠ During operation of the analyzer, the user must independently monitor the correctness of the current date and time and correct them if necessary.

2.6.3 If you need to print measurement reports on a thermal printer, prepare:

- check and, if necessary, adjust the analyzer's print settings in accordance with 2.5.7 of this manual;

- prepare the thermal printer for operation in accordance with 2.7.2 of this OM.

2.6.4 If it is necessary to transfer measurement protocols to an external device via Bluetooth prepare:

– make sure that the analyzer is located at a distance of no more than 10 m from the external device, and there is no interference with the transmission of the measurement protocol (according to 2.1.16 of this manual)

– make sure that the analyzer has the following print settings (according to 2.5.7 of this manual):

- auto print function activated;
- one copy of the protocol has been installed for one and the other version of the testimony;

– make sure that the analyzer is connected to an external device (the Bluetooth indicator is displayed in the main menu of the analyzer).

2.6.5 If you need to transfer measurement protocols to an external device via an interface cable, make sure that the analyzer has the following print settings (according to 2.5.7 of this manual):

- auto print function activated;
- There is one copy of the protocol for one and the other version of the testimony.

2.7 Preparing and operating a thermal printer

2.7.1 Installing and Removing the Battery Pack from a Thermal Printer

2.7.1.1 To install the battery pack into the thermal printer, do the following:

a) open the cover of the thermal printer's power compartment by pressing the lever on the cover;

b) insert the battery pack into the power compartment of the thermal printer; to do this, first align the contacts of the battery pack and the thermal printer, and then press the battery pack until it fits snugly into the power compartment;

c) install the thermal printer power compartment cover; to do this, first align the protrusions on the cover with the internal grooves in the thermal printer body, then press the cover until it clicks.

2.7.1.2 To remove the battery pack from the thermal printer, do the following:

a) open the cover of the thermal printer's power compartment by pressing the lever on the cover;

b) pull the battery pack up, holding the special protrusion on the battery pack (there is a recess for it in the power compartment of the thermal printer), and remove the battery pack from the power compartment;

c) install the thermal printer power compartment cover; to do this, first align the protrusions on the cover with the internal grooves in the thermal printer body, then press the cover until it clicks.

2.7.2 External inspection and preparation for operation of the thermal printer

Before starting work with the thermal printer, carry out its external inspection and preparation for work:

– check the presence and correct installation of thermal paper in the thermal printer in accordance with 3.4.4 of this manual;

– check the serviceability of the thermal printer, to do this, print the measurement protocol from the analyzer's memory in accordance with 2.13 of this manual in order to assess the print quality (prints of characters in the measurement report on paper must be easy to read, the information corresponds to the manual and the passport of the analyzers, the serial number of the analyzer is indicated correctly);

– make sure that the thermal printer is located no more than 10 meters from the analyzer;

– make sure that there is no interference with the transmission of the measurement protocol (in accordance with the instructions of 2.1.16 of this OM);

– if necessary, charge the battery pack (in accordance with 3.4.3 of this OM) or, if applicable, connect the thermal printer to the power supply in accordance with 1.2.16.1 of this OM.

2.7.3 How to use a thermal printer

2.7.3.1 Turn on the thermal printer by long (about 2-3 seconds) pressing the on / off button on the thermal printer, and the on / off button indicator light will light up.

Note – Depending on the modification of the thermal printer, after turning it on, the indicator light of the on / off button may begin to blink. The blinking frequency of the indicator light decreases (or stops blinking) to indicate that a Bluetooth connection has been established.

2.7.3.2 If after 5-10 seconds the connection between the analyzer and the thermal printer is not automatically established (the Bluetooth indicator does not appear in the main menu of the analyzer), then carry out the connection according to 2.5.8 of this manual.

 **If the thermal printer is not used for more than 5 minutes, it automatically turns off in order to save battery power, and the on / off**

button indicator light goes out. If you need to continue working with the thermal printer, turn it on again.

2.7.3.3 If markings on paper in the form of colored stripes appear in the printed measurement protocol, indicating that the paper in the thermal printer is running out, replace the paper in the thermal printer in accordance with 3.4.4 of this OM.

2.7.3.4 When finished, turn off the thermal printer by long (about 3 seconds) pressing the on / off button, and the on / off button indicator light will go out.

2.8 Operating procedure in measurement mode with data entry

2.8.1 Before taking measurements, make sure that the conditions listed in 2.1 of this OM are met, prepare the analyzer for operation in accordance with 2.6.2 of this OM, familiarize the subject with the procedure for taking measurements, as well as the rules for exhalation (in accordance with the instructions given in 2.8.13 of this RE).

2.8.2 Turn on the analyzer in accordance with 2.2 of this manual.

2.8.3 Set the measurement mode with data entry in accordance with 2.5.5.4 of this manual.

2.8.4 In the main menu, select the measurement mode by clicking on the corresponding icon .

A data entry line and a virtual keyboard for typing characters will appear on the screen (Figure 14).

Note – If the analyzer's verification results have expired, and the verification work blocking function is activated in the analyzer (in accordance with 1.4.9 of this OM), then when selecting the measurement mode with data entry, the message "Requires Verification" is displayed on the screen; For further operation, the analyzer must be submitted for verification.

2.8.5 Enter data using the virtual keyboard.

When you click on the virtual button "АБВ" (or "ABC"), the virtual keyboard layout changes from Cyrillic to Latin with numbers and back (Figure 14).

Cyrillic

Latin alphabet, numbers

Figure 14 – Data entry windows in measurement mode with data entry (input field for the name of the subject)

When you press the virtual button "DELETE" (or "DEL"), the last typed character is deleted from the input line.

2.8.6 To move to the next field, press virtual button , and if you leave the fields blank (for example, if there is no data to enter), then you can enter the missing data into the printed measurement report by hand.

Note – The entered data of fields 2, 4, 5 and 6 – option 1 of factory settings (according to 1.4.11 of this OM) or fields 3, 5 and 6 – option 2 of factory settings (according to 1.4.12 of this OM) are stored in the analyzer's memory, and when entering data before measurement, the previously entered field value is immediately displayed in the input line for these fields. Previously saved data can be corrected by deleting the previous value and entering a new one, which will also be stored until the next change.

2.8.7 You can return to the previous field by pressing the virtual button .

 *If you return to the previous field, the entered data will be lost, except for those data stored in memory.*

2.8.8 When the data in all fields has been entered, a window with the message "Insert a clean mouthpiece" will appear on the screen (Figure 15, left).


Figure 15 – Beginning of the measurement procedure

2.8.9 After making sure that the primary packaging of the mouthpiece is not damaged, open the primary packaging of the mouthpiece.

! If the original packaging is damaged, do not use the mouthpiece.

Note – When carrying out several consecutive measurements of the mass concentration of ethanol vapor in the exhaled air of one subject, it is allowed to reuse the mouthpiece.

2.8.10 **i** Install the mouthpiece. To do this, first, holding the mouthpiece with the holes down at an angle of 45°, place the locking protrusion of the mouthpiece under the mouthpiece grip in the upper part of the analyzer body (Figure 1, Figure 2), then lower the mouthpiece until it stops.

i Looking at the side panel of the analyzer (from the side of the mouthpiece pusher), visually check that the mouthpiece is installed correctly: both analyzer fittings should be completely recessed into the corresponding holes in the mouthpiece.

2.8.11 Press the virtual “Start” button on the analyzer screen to begin the measurement procedure.

When you press the “Start” button, an air sample is taken to automatically check for the absence of ethanol and other components that may affect the analyzer readings in the ambient air, in the analyzer intake system and in the mouthpiece (Figure 15, right):

– if the test result is positive (in the absence of ethanol and other components), the analyzer will go to the window ready for sampling (Figure 16, right);

Alcohol detected
will appear on the screen, in this case proceed according to 2.8. 20 of this RE.


Figure 16 – Readiness for sampling

2.8.12 When the analyzer is ready to take a sample, exhale in accordance with 2.8.13 of this OM.

If the exhalation was not done in a timely manner (within 30 seconds from the moment the window indicating readiness for sample collection appears), then to return the analyzer to the mode of readiness for sample collection, proceed according to 2.8. 21 of this RE.

To interrupt the measurement procedure in order to record the fact of refusal to take measurements (when the subject refuses to exhale), proceed in accordance with 2.10 of this OM.

2.8.13 The person being examined should blow into the inlet of the mouthpiece (Figure 2).

The exhalation should be calm and uniform (without forcing) and last about 4-5 seconds (until the end of the reference sound signal). If the subject provides the required flow of exhaled air (at least 20 l/min), a reference sound signal sounds during exhalation, and the remaining exhalation time as a percentage is shown on the exhalation duration indicator (Figure 1 7). At the end of the reference sound signal, the analyzer will

automatically take a sample of exhaled air, after which it is advisable to stop exhaling.

Automatic sampling is performed at the very end of exhalation to ensure analysis of air from deep in the lungs.


Figure 17 – Exhalation indication

2.8.14 If the exhaled air sample does not meet the required parameters for flow or volume (specified in 1.2.5 of this OM), i.e. in case of interruption of exhalation, the message “**Exhalation Interrupted**” is displayed on the screen. In this case, it is advisable to repeat the exhalation in accordance with 2.8.22 of this OM.

To interrupt the measurement procedure in order to record the fact of interruption of exhalation, act in accordance with 2.8.23 of this OM.

2.8.15 If the examined person, for physiological reasons, is unable to exhale with the required parameters for flow or volume (specified in 1.2.5 of this OM), it is possible to use manual sampling in accordance with 2.9 of this OM.

2.8.16 After taking an air sample, the mass concentration of ethanol vapor in the selected sample is measured, and the message “**Analysis in progress, please wait...**” is displayed on the screen (Figure 18, left).

After the sample has been analyzed, the result of measuring the mass concentration of ethanol (alcohol) vapor in the selected sample of exhaled air is presented on the screen (Figure 18, right).


Figure 18 – Presentation of measurement results

If the automatic printing function is activated (in accordance with 2.5.7 of this manual), then after the measurement result is displayed on the screen, printing the measurement protocol on a thermal printer or transferring the measurement protocol to an external device automatically begins.

2.8.17 If you need to print a measurement protocol (or transfer it to an external device), click on the virtual button .

Note – The measurement protocol is printed on a thermal printer in accordance with the print settings of the analyzer (in accordance with 2.5.7 of this manual).


2.8.18  Remove the used mouthpiece and dispose of it in accordance with the instructions of 4.3.2 of this manual.

2.8.19 To carry out the next measurement, press the virtual button  and repeat the procedure according to 2.8.5-2.8.18 of this manual.

To exit to the main menu, click on the virtual button .

To complete the work, turn off the analyzer by long (about two seconds) pressing the on/off button until the “**Shutdown**” message appears on the screen.

2.8.20 If, when performing 2.8.11 **Alcohol detected** ” will appear on the screen (Figure 19, right).

Ethanol vapors were detected in the ambient air, mouthpiece or intake system

To enter sample ready mode, do the following:

- return to the beginning of the measurement procedure by pressing the virtual button , in this case all data entered before the measurement will not be lost;
 - make sure that the operating rules according to 2.1.4 of this manual are observed, since ethanol vapor could enter the mouthpiece from the ambient air;
 - repeat the procedure according to 2.8.11 of this OM;
- Alcohol detected** ” appears repeatedly, remove the mouthpiece and repeat the procedure according to 2.8.9-2.8.11 of this manual.

 *As long as ethanol vapor is detected in the analyzer's mouth-piece or intake system, the analyzer will not be ready to sample.*

2.8.21 If, when performing 2.8.12 **Blow** ” window of readiness for sampling (Figure 20, left) will be replaced by the “ **No exhalation** ” window (Figure 20, right) after 30 seconds.


Figure 20 - Changing windows in the absence of exhalation

To resume the measurement procedure, press the virtual button “ **Repeat** ”, in this case all data entered before the measurement will not be lost.

2.8.22 If, when performing 2.8.13 of this OM, the exhaled air sample does not meet the required parameters for flow or volume (specified in 1.2.5 of this OM), i.e. if exhalation is interrupted, the “ **Exhalation Interrupted** ” window will be displayed on the screen (Figure 21, right).


Figure 21 – Breathing test failure

It is advisable to stop exhaling. To resume the measurement procedure, click on the “ **Repeat** ” virtual button ; after the analyzer enters the mode of readiness for sampling (Figure 1-6 , right), you can repeat the exhalation in accordance with 2.8.13 of this manual.

2.8.23 In case of interruption of exhalation, the analyzer provides the ability to interrupt the measurement procedure, while recording the fact(s) of interruption of exhalation (if necessary).

To do this, when the “**Exhalation Interrupted**” window is displayed on the screen (Figure 21, right), press the virtual button “**Exhalation Interrupted**” - a corresponding message will appear on the screen (Figure 22).


Figure 22 – Result “Exhalation Interrupted”

the message “**Exhalation Interrupted**” as a measurement result, which indicates that an attempt was made to exhale, which did not result in automatic sampling with subsequent analysis, but in the printed (or transmitted to external devices) measurement protocol instead of the measurement result, the message “**Exhalation Interrupted**” will be displayed (there will be no line with the sampling mode).

To continue work, follow instructions 2.8.17-2.8.19 of this OM.

2.9 Manual sampling mode

2.9.1 If the subject is unable to exhale, at which automatic sampling is triggered, it is advisable to use the manual sampling mode.

To do this, during the exhalation of the subject, press the virtual button “**MANUAL SAMPLING**” (Figure 17), and the analyzer will take a manual sample. Perform manual sampling at the very end of the exhalation that the subject is capable of.

2.9.2 The successful use of exhaled air as a diagnostic medium for determining the alcohol content in a person's blood is made possible by the fact that the amount of alcohol evaporating into the exhaled air depends on its content in the venous blood. This dependence is a consequence of Henry's law, which states that the pressure of a volatile (gaseous) component

in a solution and, therefore, the vapor content of this component above the solution is directly proportional to its content in the solution.

The closest contact with the blood is pulmonary air, which is close to alveolar air, the analysis of which becomes possible only with deep exhalation (volume of at least 1.2 liters). The air in the oral cavity, as well as in the upper part of the respiratory tract, is not in close contact with venous blood.

That is why, to measure the mass concentration of ethanol vapor in exhaled air, it is recommended to use automatic sampling, which ensures compliance with the parameters of the respiratory test (exhaled air flow rate and volume).

In manual sampling mode, the exhalation parameters specified in 1.2.5 of this OM are not provided.

When manually sampling, the result of measuring the mass concentration of ethanol in the selected sample of exhaled air may be slightly underestimated⁷ compared to the result obtained when automatic sampling is triggered, due to the fact that during manual sampling, the sample contains a larger proportion of air from the upper respiratory tract. The air from the upper respiratory tract is less saturated with ethanol vapor than the pulmonary air, which is sampled during automatic sampling.

2.10 Refusal to measure

2.10.1 **Blow**” window of readiness for sampling will be replaced by the “**No exhalation**” window (Figure 23, left).

2.10.2 If the subject refuses to take measurements (the subject refuses to exhale), the analyzer provides the ability to interrupt the measurement procedure, while recording the fact of refusal to take measurements (if necessary). To do this, click on the virtual button “**Refusal**” - a message about refusal of measurements will appear on the screen (Figure 23, right).

⁷Provided there is no ethanol present on the oral mucosa.


Figure 23 – Result of “Refusal of measurements”

Refusal as a measurement result, which indicates that there was no exhalation of the subject, and in the printed (or transmitted to external devices) measurement protocol, instead of the measurement result, the inscription “Refusal of the test” will be displayed. (in this case, there will be no line with the sampling mode).

To continue work, follow instructions 2.8.17-2.8.19 of this OM.

2.11 How to operate in measurement mode without data entry

2.11.1 Before taking measurements, make sure that the conditions listed in 2.1 of this OM are met, prepare the analyzer for operation in accordance with 2.6.2 of this OM, familiarize the subject with the procedure for taking measurements, as well as the rules for exhalation (in accordance with the instructions given in 2.8.1 3 of this OM).

2.11.2 Turn on the analyzer in accordance with 2.2 of this manual.

2.11.3 Set the measurement mode without data entry in accordance with 2.5.5.4 of this OM.

In the main menu, select the measurement mode by clicking on the corresponding icon . A window with the message “Insert a clean mouthpiece” will appear on the screen (Figure 24, left).


Figure 24 – Starting the measurement procedure in measurement mode without entering data

Note – If the analyzer’s verification results have expired, and the verification blocking function is activated in the analyzer (in accordance with Notes 1.4.9 of this OM), then when selecting the measurement mode without entering data, the message “Requires Verification” is displayed on the screen . , for further operation the analyzer must be submitted for periodic verification.

2.11.4 Carry out the measurement in accordance with 2.8.9-2.8.19 of this OM.

2.11.5 In the fields of the printed measurement protocol intended for filling out, it is possible to enter data by hand.

2.12 Operating procedure in screening mode

2.12.1 Before testing, make sure that the conditions listed in 2.1 of this OM are met, prepare the analyzer for operation in accordance with 2.6.2 of this OM, familiarize the subject with the testing procedure, as well as the rules for exhalation (in accordance with the instructions given in 2.12.5 of this RE).

2.12.2 Turn on the analyzer in accordance with 2.2 of this manual.

2.12.3 Select screening mode .

A message indicating that you are ready for the test will appear on the touch screen (Figure 25, left).


Figure 25 – Screening mode

Note – If the analyzer's verification results have expired, and the verification work blocking function is activated in the analyzer (in accordance with notes 1.4.9 of this OM), then when selecting the screening mode, the message "Requires Verification" is displayed on the screen, for further operation, the analyzer must be submitted for verification.

2.12.4  Install the funnel mouthpiece as follows: holding the funnel mouthpiece with the holes down at an angle of 45°, place the fixing protrusion of the funnel mouthpiece under the mouthpiece grip in the upper part of the analyzer body (Figure 1, Figure 3), then lower the funnel mouthpiece until it stops.

 Looking at the side panel of the analyzer (from the mouthpiece pusher side), visually check that the funnel mouthpiece is installed correctly: both analyzer fittings should be completely recessed into the corresponding holes in the funnel mouthpiece.

2.12.5 To conduct testing, the subject must blow into the funnel mouthpiece strongly and continuously (but without forcing) for about two seconds. The funnel mouthpiece should be located directly opposite the subject's mouth, almost perpendicular to the surface of the face, at a distance of 1.5-2 cm.

Note – For better operation of automatic sampling, it is recommended to direct the flow of exhaled air into the elongated part of the funnel mouthpiece (opposite the hole in the funnel mouthpiece to the pressure sensor - marked with an arrow in Figure 3).

2.12.6 If the subject provides the required flow of exhaled air, a reference sound signal sounds during exhalation, and the remaining exhalation time as a percentage is shown on the exhalation duration indicator (Figure 25, right). It is advisable to stop exhalation at the end of the reference sound signal. The analyzer will automatically take a sample of exhaled air from the funnel mouthpiece.

 **Do not remove the analyzer from the subject's mouth until the sample is taken. Otherwise, the ambient air will be sampled rather than the exhaled air.**

2.12.7 If the examined person, for physiological reasons, is unable to exhale with the required flow or volume parameters, it is advisable to use manual sampling in accordance with 2.9 of this OM.

2.12.8 Depending on the test result, one of two messages appears on the analyzer screen (Figure 26).

 **The screening mode can only be used to assess the presence of ethanol in exhaled air, and the analyzer error may exceed the limits of the permissible analyzer error specified in 1.2.1 of this OM.**


Figure 26 – Test results in screening mode

The message "detected" or "not detected" (depending on the test result) is stored in the analyzer's memory as a measurement result.

2.12.9 After the test result is displayed on the screen, the result is automatically transferred to an external device.

2.12.10 If alcohol is detected in the exhaled air, it is recommended to carry out the measurement in one of the measurement modes (in accordance with 2.8 or 2.11 of this OM).

2.12.11 To conduct the next test, click on the virtual button  and repeat the procedure according to 2.12.5-2.12.10 of this manual.

 *When using the analyzer for repeated measurements using a reusable funnel mouthpiece, especially at low temperatures, it is necessary to periodically wipe the inner surface of the funnel mouthpiece (in accordance with 3.4.6.3 of this manual), preventing the appearance of condensation. Condensation entering the intake system can lead to failure of the electrochemical sensor.*

To exit to the main menu, press the virtual button .

To complete the work, turn off the analyzer by long (about two seconds) pressing the on/off button until the “Shutdown” message appears on the screen and remove the funnel mouthpiece.

2.13 View and print (transfer to an external device) measurement reports stored in the analyzer memory

2.13.1 The main menu item “View test memory” allows you to display all measurement results stored in memory.

2.13.2 View test memory” icon in the main menu of the analyzer.

Record” window will appear on the touch screen with the protocol of the last measurement or the last viewed protocol (Figure 27).


Figure 27 – Viewing test memory

2.13.3 If measurements were made in measurement mode with data entry, the measurement results are stored in memory together with automatically saved data (in accordance with 1.4.7 of this manual) and with data entered into the fields intended for filling in (Figure 27, left).

2.13.4 If measurements were made in the measurement mode without entering data, only the measurement results are stored in memory together with automatically saved data (according to 1.4.7 of this manual) (Figure 27, right).

2.13.5 Iterating through records one at a time is done by pressing the virtual buttons  and .

2.13.6 Iterating through 10 records is done by briefly pressing the service button (10 records back) or briefly pressing the on/off button (10 records forward).

Note – Viewing records in the analyzer’s memory (after one or 10 tests) is carried out according to the numbers of measurement protocols, without taking into account the date and time of their implementation. It should be taken into account that after the analyzer has carried out more than 10,000 measurements (tests), subsequent protocols are recorded in the analyzer’s memory instead of those saved, starting with number 00001. For example, it may turn out that after the measurement protocol with number 00001 with the date of December 15, 2021, the next In the analyzer memory there will be a measurement protocol with number 00002 with a date of October 10, 2020 (protocol data with number 00002 will be replaced with new ones only after another measurement).

2.13.7 To quickly jump to the first, to the last record, as well as to the protocol record under the required number, do the following:

– while in the “Record” window (Figure 27), press and hold the service button and briefly press the on/off button, release both buttons; a window for switching to the required protocol will appear on the screen (Figure 28);

– to go to the first or last record, click on the virtual button “Start” or “End”, respectively;

Note – After the analyzer has carried out more than 10,000 measurements, pressing the “Start” button switches to protocol number 00001, and pressing the “End” button switches to the last measurement protocol in memory.

– to quickly go to a measurement protocol with a specific number, enter the required test number using the “+” or “-” buttons (the most significant digit of the number is indicated on the left), and click on the “Go” virtual button. If the test with the dialed number does not exist, the transition is made to the last measurement protocol in memory.


Figure 28 – Quick transition to the required protocol

2.13.8 In the test memory viewing menu, depending on the connection status of the analyzer with a thermal printer (or an external device via Bluetooth), the title of the “Recording” window displays the words “Connected” or “Disconnected”:

Connected status – the analyzer is ready to print or transfer the measurement protocol to an external device via Bluetooth;

Disabled status – the analyzer is not ready to print or transfer the measurement protocol to an external device via Bluetooth.

Note – If in the menu for viewing the test memory in the title of the “Recording” window the message “Disconnected” is displayed due to the lack of connection with a thermal printer, then it is impossible to print a measurement protocol. In this case, it is necessary to connect the analyzer to the thermal printer in accordance with 2.5.8 of this manual.

2.13.9 When you press the virtual button, the measurement protocol displayed on the screen is printed on a thermal printer or transferred to an external device (via Bluetooth or interface cable), provided there is a connection.

Note – The virtual “Print” button is absent in the test memory viewing window when displaying the results of tests performed in screening mode.

2.14 Contents of the printed measurement report

2.14.1 In the analyzer measurement protocols, information is printed according to Table 9 (an example is given with a set of fields for data entry in accordance with 1.4.11 of this OM).

Table 9 – Contents of the measurement protocol

Inscription in the protocol	Contents of the protocol
AJKOTEKTOP PRO-100 touch-M plus	Analyzer type designation
Device Number:	Analyzer serial number
Adjustment date: DD/MM/YYYY	Date of last adjustment of analyzer readings (day/month/year)
Verification date: DD/MM/YYYY	Date of the last analyzer verification (day/month/year)
Test NO.:	Measurement number (according to the internal numbering of the analyzer)
Date: DD/MM/YYYY	Date of measurement (day/month/year)
Time: HH:MM	Measurement time (hour:minutes)
Mode: Automatic	Air sampling mode ¹⁾
Result: X.XXX mg/l	The result of measuring the mass concentration of ethanol vapor in a selected sample of exhaled air: numerical value and designation of the unit of measurement “mg/l” ²⁾
Subject's name: ⁶⁾	Last name and initials of the person being examined ³⁾
Survey Location: ⁶⁾	Measuring location data ³⁾
State Vehicle Number: ⁶⁾	State number of the vehicle ³⁾
Badge: ⁶⁾	Inspector badge number ³⁾
Inspector: ⁶⁾	Inspector's surname and initials ³⁾
Traffic police department: ⁶⁾	Traffic police department number ³⁾
Coordinates:	Coordinates of the measurement location ⁴⁾
Signature Survey: ⁷⁾	Signature of the person being examined ⁵⁾
Signature: ⁷⁾	Inspector's signature ⁵⁾

End of table 9

- ¹⁾ In the manual air sampling mode, the message "Mode: MANUAL SAMPLING" is displayed.
- ²⁾ If it is recorded that the subject refuses to take the measurement, the message "Refusal of the test" is displayed; in the event of a recorded fact of insufficient flow and exhalation volume, the message "Exhalation Interrupted" is displayed. In this case, information about the air sampling mode is not displayed in the protocol.
- ³⁾ Data is entered from the analyzer's virtual keyboard (touch screen) before measurement or entered by hand into a printed measurement protocol.
- ⁴⁾ The data is printed in the measurement report when the coordinates of the measurement location are established and the function of entering coordinates into the measurement report is activated; metrological characteristics for the coordinates of the measurement location are not standardized.
- ⁵⁾ The data is entered by hand into the printed measurement protocol.
- ⁶⁾ The set of fields for data entry may differ from those indicated in the table (the number of fields is set from 0 to 10, the names of the fields can be changed). The names of the fields for entering measurement protocol data are indicated in the analyzer passport.
- ⁷⁾ The fields are printed in the measurement report when the function of entering fields for signature in the measurement report is activated.

2.14.2 When transmitting a measurement result to an external device via Bluetooth or an interface cable, the transmitted protocol contains information according to Table 9 (except for signature fields - these fields are not transmitted).

2.15 Turning off the analyzer

2.15.1 The analyzer is turned off by long (about three seconds) pressing the on/off button until the "Shutdown" message appears on the screen.

2.15.2 In order to save battery power, the analyzer has an automatic shutdown function. The user can independently set a convenient time for automatic shutdown (according to 2.5.6 of this manual).

2.16 Turning off the thermal printer

2.16.1 on / off button on the thermal printer for a long time (about three seconds), and the on / off button indicator light goes out.

2.16.2 If the thermal printer is not used for more than 5 minutes, it automatically turns off to save battery power.

2.17 Possible malfunctions of the analyzer and thermal printer

Possible malfunctions of the analyzer and thermal printer, as well as methods for eliminating them, are given in Table 10.

Table 10 - Possible faults

No.	Possible malfunction	Probable cause	Way elimination
1	The analyzer does not turn on	1. No power supply	Install/charge batteries in accordance with 3.4.2 of this OM or connect the analyzer to the power supply in accordance with 1.2.15 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
2	The real time clock has stopped (the current date and time are not set)	1. Lithium battery undervoltage 2. The analyzer is faulty	The analyzer must be sent to a service center*
3	The analyzer does not automatically take a sample	1. Insufficient flow or volume of exhaled air	Repeat exhalation in accordance with 2.8.13 or 2.12.5 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
4	There is no sound signaling of the analyzer operation	The analyzer is faulty	The analyzer must be sent to a service center *

Continuation of table 10

5	The measurement report prints out faint text or no text at all	1. Thermal printer battery pack voltage is low	Charge the battery pack in accordance with 3.4.3 of this OM or, if applicable, connect the AC adapter charger to the thermal printer
		2. Poor paper pressure with the thermal printer's rubber roller	Make sure the thermal printer's pressure roller is clean; clean it if necessary; Close the thermal printer cover tightly. If the pressure roller is damaged, replace it**
		3. The thermal paper roll is installed on the wrong side, or the wrong paper is being used.	Make sure that the correct paper is used and check that it is installed correctly according to 3.4.4 of this manual
		4. The thermal printer is faulty	The thermal printer must be sent to a service center *
6	The measurement protocol is not printed after the measurement result is displayed on the analyzer screen	1. Incorrect print settings	Set the required number of copies according to 2.5.7 of this OM
		2. There is no connection with the thermal printer via Bluetooth	Establish a connection between the analyzer and the thermal printer in accordance with 2.5.8 of this manual
		3. Presence of interference with the transmission of the measurement protocol via Bluetooth (specified in 2.1.16 of this OM)	Repeat the transfer by clicking on the virtual button 

Continuation of table 10

		4. The analyzer and/or thermal printer is faulty	The analyzer together with the thermal printer must be sent to the service center *
7	The measurement protocol is not transmitted to an external device via Bluetooth after the measurement result is displayed on the analyzer screen	1. Incorrect print settings	Activate the automatic printing function and install one copy in accordance with 2.5.7 of this OM
		2. Presence of interference with the transmission of the measurement protocol via Bluetooth (specified in 2.1.16 of this OM)	Repeat the transfer by clicking on the virtual button 
		Bluetooth connection	Make sure the Bluetooth connection is established (the Bluetooth indicator appears in the main menu). Check the operation of the external device
8	The measurement protocol is not transmitted to an external device via the interface cable after the measurement result is displayed on the analyzer screen	1. Incorrect print settings	Activate the automatic printing function and install one copy in accordance with 2.5.7 of this OM
		2. An interface cable model is used that does not support data transfer to an external device.	Make sure you are using the correct model of interface cable for the individual analyzer functionality listed in its data sheet.
		3. The analyzer and/or interface cable is faulty	The analyzer along with the interface cable must be sent to the service center *

End of table 10

9	When turning on the analyzer, the message "Temperature very high (low)" appears and the analyzer turns off	1. Temperature of the measuring board is more than +50 °C (less than -5 °C)	Maintain the analyzer under operating conditions (specified in 1.1.6 of this manual) and turn it on again
		2. The analyzer is faulty	The analyzer must be sent to a service center *
10	After pressing the "Start" button, the analyzer does not enter the window ready for sample collection; the message "Alcohol detected" appears on the screen.	1. Presence of ethanol in the mouthpiece or intake system of the analyzer	Act in accordance with 2.8.20 of this OM
		2. The analyzer is faulty	The analyzer must be sent to a service center *
eleven	The analyzer does not respond correctly or does not respond to pressing virtual buttons on the touch screen	1. Touch screen calibration reset	Calibrate the touch screen in accordance with 3.4.5 of this manual
		2. The analyzer is faulty	The analyzer must be sent to a service center*

* During the warranty period, the analyzer should be sent to ALKOTEKTOR LLC, during the post-warranty period - to ALKOTEKTOR LLC or to other organizations that have the appropriate licenses and permits (service centers). Information about service centers is given in 6.3 of this OM.

** Pressure roller for thermal printer is available upon request.

3 MAINTENANCE

3.1 Maintenance of the analyzer is carried out in order to maintain the analyzer in good condition and ensure its constant readiness for use.

3.2 When carrying out maintenance of the analyzer, in order to ensure safety, including environmental safety, the requirements of regulatory documents on occupational health and safety, including those contained in

the operational documentation of the tools, equipment and devices used during maintenance, must be observed.

3.3 Analyzer maintenance includes:

- on-site maintenance;
- periodic maintenance.

3.4 On-site maintenance

3.4.1 On-site maintenance is carried out by the user as necessary, including based on the results of checking the analyzer before starting work, which is carried out during an external inspection and preparation of the analyzer and thermal printer for operation in accordance with 2.6.2 and 2.7.2 of this OM, and includes:

- charging/replacing batteries at low voltage levels;
- charging the battery pack at low voltage;
- refilling the thermal printer with thermal paper;
- calibration of the touch screen - as necessary (in case the analyzer does not respond correctly or does not respond to pressing virtual buttons/icons on the touch screen);
- cleaning and disinfection of the analyzer as it gets dirty (frequency is set by the user).

3.4.2 Charge /replace analyzer batteries

3.4.2.1 To check the supply voltage level, go to the main menu of the analyzer.

When the battery voltage level is low, only a small part of the power indicator in the main menu remains shaded; When the battery voltage level drops below the permissible level, the analyzer will turn off. In these cases, it is necessary to charge/replace the batteries.

Note – The analyzer shutdown in the event of low power can occur at all stages of its operation, so monitor the power indicator in the main menu and charge/replace batteries in advance.

3.4.2.2 To replace/charge batteries, follow these steps:

- a) turn off the analyzer;
- b) disconnect the power compartment cover by sliding it in the direction of the arrow on the cover;
- c) carefully pull the protruding part of the auxiliary tape up and remove all 4 batteries from the compartment;
- d) install new batteries/charged accumulators into the power compartment, observing the polarity; at the same time, make sure that the tape

passes under the batteries, and its remaining part protrudes between the outermost battery and the wall of the power compartment;

e) install the power compartment cover, aligning the tabs on the cover with the internal grooves on the analyzer body, and push the power compartment cover until it stops.

 **Install batteries of the same type and manufacturer in the power compartment. When replacing batteries, replace all 4 batteries at the same time.**

3.4.3 Charging the battery pack (for thermal printer)

3.4.3.1 If your thermal printer has a power connector to charge the battery pack, follow these instructions:

a) install the battery pack into the power compartment of the thermal printer (in accordance with 2.7.1.1 of this manual);

b) connect the network or on-board power adapter of the charger to the power connector of the thermal printer;

c) connect the charger power adapter to the power supply; the charge must be continued until the indicator light on the on / off button of the thermal printer goes out ;

d) disconnect the charger power adapter from the power supply, then from the thermal printer.

3.4.3.2 If you have a charger, follow these steps to charge the battery pack:

a) remove the battery pack from the power compartment of the thermal printer (in accordance with 2.7.1.2 of this manual);

b) install the battery pack into the charger: first align the contacts of the battery pack with the contacts of the charger, then lower the battery pack down until it stops;

c) connect the charger power adapter to the charger;

d) connect the charger power adapter to the power supply, the charger indicator will turn red; the charge must be continued until the indicator color changes to green;

e) install the battery pack back into the power compartment of the thermal printer.

3.4.4 Refilling the thermal printer with thermal paper

3.4.4.1 To load a thermal paper roll into the thermal printer, follow these steps:

a) pull the tabs on the thermal printer cover towards you;

b) insert a new roll of thermal paper, positioning the roll as shown in Figure 29;

c) close the thermal printer cover by pressing it until it clicks.


Figure 29 – Appearance of the paper compartment cover of a thermal printer

Note – The thermal side of the paper is easy to determine: if you quickly scratch the stylus on the thermal side, a black mark will remain.

3.4.5 Touch screen calibration

3.4.5.1 To enter the analyzer touch screen calibration mode, follow the instructions:

a) press the analyzer on/off button and immediately, without releasing it, briefly press the service button (within 0.5 seconds);

b) a white screen will appear on the analyzer screen with a warning about entering the touch screen calibration mode;

c) touch any point on the screen with the stylus - the message “**CALIBRATION 3 touches**” will be displayed in the center of the screen left”;

d) touch the points with the stylus, which will be indicated in turn on the analyzer screen;

Note – Touching the analyzer’s touch screen with the stylus should be done gently, but with sufficient force (for example, as when writing with a pen on paper).

e) after the screen calibration is carried out, the message “**Touch screen**” will be displayed on the screen to continue”;

f) touch any point on the screen with the stylus - the analyzer will exit the touch screen calibration mode, and the main menu of the analyzer will appear on the screen.

3.4.6 Cleaning and disinfection

3.4.6.1 Clean the analyzer housing, thermal printer and power adapters with a slightly damp cloth.

Do not use abrasive chemicals for cleaning as this may damage the housing, touch screen and/or electrochemical sensor of the analyzer.

When cleaning, make sure there is no dirt or dust in the inlet ports, analyzer intake system. If present, contamination should be removed.

3.4.6.2 Clean the rubber pressure roller of the thermal printer with a dry, lint-free cloth.

3.4.6.3 To remove condensate from the funnel mouthpiece of the analyzer during operation, use a dry, lint-free cloth (not impregnated with alcohol-containing or other chemical cleaning agents).

3.4.6.4 The analyzer must be periodically sanitized according to MU-287-113 in accordance with the standards established V institution where _ the analyzer is operated (frequency is set according to the degree of contamination, similar to how For medical tools , not having contact with mucous shells body person , in compliance With profile medical institutions).

Disinfection should be carried out using the chemical method according to MU-287-113. For disinfection, the method of wiping twice with a calico or gauze napkin soaked in a disinfectant solution and wrung out to prevent the disinfectant solution from getting inside the analyzer can be used. For disinfection, it is allowed to use chemicals recommended for plastic products, except for alcohol-containing disinfectants , for example, a 3% solution of hydrogen peroxide with the addition of 0.5% Lotus detergent .

 After disinfection, the funnel mouthpiece should be thoroughly dried naturally. method to prevent disinfectant from entering the analyzer.

 *M mouthpieces disinfection Not are exposed , and are disposed of after application in accordance with 4.3.2 of this OM.*

3.5 Periodic Maintenance

3.5.1 Periodic maintenance of the analyzer is performed at the service center during the entire period of operation and includes:

– technical condition monitoring – recommended before verification;

– checking the analyzer readings ⁸ – it is recommended before verification (or more often in accordance with the frequency set by the user);
– adjustment of analyzer readings ⁸ – in case of negative results of checking the readings or during verification as necessary.

Notes:

1 When calibrating the analyzer, the readings are adjusted if necessary in accordance with the instructions in paragraph 8.3.3 of the document MP- 242-2554-2023 "Analyzers of ethanol vapors in exhaled air ALKOTEKTOR PRO-100 touch - M plus . Verification methodology."

2 Checking and adjusting readings can be performed by authorized technical service specialists of the organization operating the analyzer, if they have the maintenance tools specified in Appendix D (except for a multimeter).

3.5.2 A description of the analyzer maintenance operations is given in the document "Analyzers of ethanol vapor in exhaled air "ALCOTEKTOR PRO-100 touch - M" with accessories. Maintenance instructions. Engineer's menu", which is available to service centers upon additional request.

Note – For organizations that have their own maintenance service and maintenance facilities, upon additional request, an extract from the document "Analyzers of ethanol vapors in exhaled air "ALCOTEKTOR PRO-100 touch -M" with accessories can be provided . Maintenance instructions. Engineer's menu", which contains all the necessary information for checking and adjusting readings.

3.5.3 The maintenance items required to perform periodic maintenance of the analyzer at service centers are listed in Appendix E.

3.5.4 It is recommended to mark records of performed maintenance and repairs in the analyzer passport in the maintenance accounting table, for example, with the entry "Checking readings", "Adjusting readings", etc.

Note - Entries in the passport about work performed are made only by service center specialists and verifiers ; If necessary, the user is recommended to keep a separate maintenance log for the analyzer.

⁸Checking and adjusting the analyzer readings can be performed as an independent maintenance operation, or as part of the procedure for monitoring the technical condition of the analyzer.

3.5.5  After correcting the analyzer readings, it is necessary to calibrate the analyzer.

4 TRANSPORTATION, STORAGE AND DISPOSAL

4.1 Transportation

4.1.1 The analyzer in transport packaging is transported by all types of transport in covered vehicles in accordance with the rules for the transportation of goods in force for each type of transport.

The placement and fastening of packaged analyzers in vehicles must ensure their stable position and exclude the possibility of impacts against each other, as well as against the walls of vehicles.

4.1.2 Transportation conditions:

- air temperature: from minus 50 °C to plus 50 °C;
- upper value of relative air humidity: 100% at plus 25 °C.

4.1.3 At a minimum, it is recommended that handling signs corresponding to the following meanings be applied to transport packaging: “Fragile, be careful!”, “Keep away from moisture.”

4.1.4 After transportation in subzero temperatures, it is recommended to keep the analyzer in operating conditions for at least an hour before turning it on. If the “Low temperature” message appears on the analyzer screen when turned on, the holding time should be increased to two hours.

4.2 Storage

4.2.1 Analyzers should be stored in closed heated rooms under the following conditions :

- air temperature: from plus 5 °C to plus 40 °C;
- air humidity, no more than: 80% at plus 25 °C.

4.2.2 It is not allowed to store analyzers in places with increased dust and gas contamination, as well as in places where alcohol-containing substances are stored in open containers or where surfaces (equipment) are treated with alcohol-containing solutions.

4.2.3 During long-term storage, remove the batteries from the analyzer and the battery pack from the thermal printer to prevent damage to the analyzer and thermal printer due to the risk of electrolyte leakage.

4.2.4 The battery pack must be stored in a charged state (the ideal state of charge is 40 %). If a battery pack is left to be stored discharged for a period of time of 3 months or longer, an irreversible drop in its capacity

may occur, and corrosion of the cells may also occur.

4.2.5  The shelf life of mouthpieces is unlimited.

4.2.6 Analyzers should be stored in packaging and placed on shelves in no more than three rows.

4.2.7 Temporary anti-corrosion protection is not required.

4.2.8 When putting the analyzer into operation after long-term storage (more than one calibration interval), periodic verification is carried out.

4.3 Disposal

4.3.1 To prevent potentially harmful effects on the environment and human health, the analyzer must not be disposed of with municipal solid waste, but should be disposed of in accordance with the rules provided for in SanPiN 2.1.3684-21 for class G.

 **Electrical and electronic devices must be disposed of through designated facilities designated by local authorities, but not with household waste.**

4.3.2  In medical institutions, used mouthpieces should be disposed of in accordance with the rules provided for in SanPiN 2.1.3684-21 for class B.

5 CHECKING THE ANALYZER

5.1 In accordance with the Federal Law of June 26, 2008 No. 102-FZ "On ensuring the uniformity of measurements" (Article 13) "Measuring instruments intended for use in the field of state regulation of ensuring the uniformity of measurements are subject to initial verification before commissioning, as well as after repair, and during operation - periodic verification. Legal entities and individual entrepreneurs using measuring instruments in the field of state regulation to ensure the uniformity of measurements are required to promptly submit these measuring instruments for verification. Verification of measuring instruments is carried out by legal entities and individual entrepreneurs accredited in accordance with the legislation of the Russian Federation on accreditation in the national accreditation system for verification of measuring instruments."

5.2 Verification of the analyzer is carried out according to the document MP-242-2554-2023 "Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M plus. Verification method".

5.3 The interval between verifications is 1 year.

5.4 The verification results are confirmed by information on the verification results entered into the Federal Information Fund for Ensuring the Uniformity of Measurements (FSIS "ARSHIN").

5.5 At the user's request or if there are requirements in regulatory documents governing the use of analyzers, when entering the verification date into the analyzer's memory, the function of blocking the analyzer's operation after the verification results have expired can be activated.

When activating the function of blocking the operation of the analyzer by verification date, it is necessary to make a corresponding mark in the maintenance accounting table of the analyzer passport.

6 MANUFACTURER'S WARRANTIES AND REPAIRS

6.1 manufacturer's warranty

ALKOTEKTOR LLC, being an authorized representative of the manufacturer, accepts claims, complaints and reviews regarding the quality of the analyzer, and also bears warranty obligations in accordance with the provisions set out in the relevant section of the analyzer passport.

6.2 Warranty repair

Warranty repairs of analyzers are carried out at the service center of ALKOTEKTOR LLC at the address: 199178, St. Petersburg, emb. Smolenskiy River, 5-7, letter A, tel. (812) 320-22-97, email, address: www.alcote.ru, email, mail: info@alcote.ru

6.3 Post-warranty maintenance and repair

Post-warranty maintenance and repair of analyzers is carried out at ALKOTEKTOR LLC, as well as at regional service centers.

The current list of service centers with which ALKOTEKTOR LLC has entered into a cooperation agreement is posted on the website www.alcote.ru in the "Service in the Regions" section.

To provide services for the maintenance and repair of analyzers, the organization must have all the necessary permits in accordance with the current legislation of the Russian Federation.

When contacting these organizations to provide maintenance or repair services, ask about the availability of permits.

APPENDIX A

Regulations

1 IEC 60601-1:2020 "Medical electrical equipment . Part 1. General requirements for basic safety and essential performance" (GOST R IEC 60601-1-2022 "Products medical electric . Part 1. General requirements security With taking into account main functional characteristics ").

3 IEC 62366-1:2015 "Medical devices – Part 1: Application of usability engineering to medical devices" (GOST R IEC 62366-1-2021 " Products medical _ Part 1. Design medical products With taking into account operational suitability ").

3 IEC 62304:2006 "Medical device software – Software life cycle processes" (GOST IEC 62304-2022 " Products medical _ Software provision . Processes life cycle ").

4 ISO 14971:2019 "Medical devices – Application of risk management to medical devices" (GOST ISO 14971-2021 " Products medical _ Application of risk management to medical devices").

5 GOST IEC 60950-1-2014 "Information technology equipment. Safety requirements. Part 1. General requirements."

6 R 50.2.077-2014 "Recommendations for metrology. State system for ensuring the uniformity of measurements. Testing of measuring instruments for the purpose of type approval. Checking software security."

7 MU-287-113 Guidelines for disinfection, pre-sterilization cleaning and sterilization of medical products.

8 Sanitary rules and regulations SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of territories of urban and rural settlements, water bodies, drinking water and drinking water supply to the population, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and carrying out sanitary and anti-epidemic (preventive) measures."

9 Federal Law of June 26, 2008 No. 102- FZ " On Ensuring the Uniformity of Measurements."

Note – When using this OM, it is advisable to check the validity of the reference documents in the public information system. If the reference document is replaced (changed), then when using this operating manual you should be guided by the replacement (changed) document. If the reference document is canceled without replacement, then the provision in which a reference to it is given applies to the part that does not affect this reference.

APPENDIX B

Electromagnetic Compatibility Information

Medical electrical equipment requires special measures to ensure electromagnetic compatibility and must be installed and commissioned in accordance with the EMC-related information provided in the operating documentation.

The use of mobile radio frequency communications may affect medical electrical products.

The use of accessories, transducers and cables other than those listed, other than transducers and cables supplied by the product or system manufacturer as replacement parts for internal parts, may result in increased electromagnetic emissions or reduced immunity of the product or system.

Electromagnetic compatibility tables

Table B.1 - Guidance and manufacturer's declaration - electromagnetic emissions

Guidance and manufacturer's declaration - electromagnetic emissions		
The system is intended for use in the electromagnetic environment specified below. The purchaser or user of the system should ensure that it is used in the specified electromagnetic environment.		
Electromagnetic Emission Test	Correspondence	Electromagnetic environment - instructions
The group to which ME belongs according to CISPR 11 (GOST R 51318.11-2006)	Group 1	The system uses RF energy only to perform internal functions. RF emissions are low and are not likely to cause interference with nearby electronic equipment
The class to which ME belongs according to CISPR 11 (GOST R 51318.11-2006)	Class B	The system is suitable for use in any location, including residential buildings and buildings directly connected to the electrical distribution network supplying residential buildings
Harmonic components of current according to IEC 61000-3-2 (GOST 30804.3.2-2013)	Class A	
Voltage fluctuations and flicker according to IEC 61000-3-3 (GOST 30804.3.3-2013)	Compliant	

Table B.2 - Guidance and manufacturer's declaration - interference immunity

Guidance and Manufacturer's Declaration - Immunity			
The system is designed for use in the electromagnetic environment specified below. The purchaser or user of the system should ensure that it is used in the specified electromagnetic environment.			
Immunity test	Test level	Compliance level	Electromagnetic environment - instructions
Electrostatic discharges (ESD) according to IEC 61000-4-2 (GOST 30804.4.2-2013)	±6 kV - contact discharge	Compliant	The floor in the room is made of wood, concrete or ceramic tiles. For floors covered with synthetic material, the relative humidity is at least 30%
	±8 kV - air discharge	Compliant	
Nanosecond impulse noise according to IEC 61000-4-4 (GOST 30804.4.4-2013)	±2 kV - for power lines	Compliant	Electrical power quality on the grid according to typical conditions in a commercial or hospital environment
	±1 kV - for input/output lines		
High energy microsecond pulse noise according to IEC 61000-4-5 (GOST R 51317.4.5-99)	±1 kV when applying interference in a wire-to-wire manner	Compliant	Electrical power quality on the grid according to typical conditions in a commercial or hospital environment
	±2 kV when noise is applied via the wire-to-ground circuit		
Voltage dips, short-term interruptions and voltage changes in the input power lines according to IEC 61000-4-11 (GOST 30804.4.11-2013)	<5% U _n (voltage dip >95% U _n) for 0.5 period	Compliant	Electrical power quality on the grid according to typical conditions in a commercial or hospital environment
	40% U _n (voltage dip 60% U _n) for 5 periods		
	70% U _n (voltage dip 30% U _n) for 25 periods		
	<5% U _n (voltage dip >95% U _n) for 5 s		
Magnetic field of industrial frequency (50/60 Hz) according	3 A/ m	Compliant	Power frequency magnetic field levels should be consistent with typical

to IEC 61000-4-8 (GOST R 50648-94)			commercial or hospital environments
NOTE: U _n - voltage level of the electrical network before the test stimulus is applied.			

Table B.3 - Guidance and manufacturer's declaration - interference immunity

Guidance and Manufacturer's Declaration - Immunity NOT INTENDED TO SUPPORT VITAL FUNCTIONS			
The system is intended for use in the electromagnetic environment specified below. The purchaser or user of the system should ensure that it is used in the specified electromagnetic environment.			
Equipment stability testing	Test level	Compliance level	Electromagnetic environment - instructions
			Portable and mobile radio frequency equipment, incl. cables should not be used near the device closer than the recommended distance, calculated by the formula according to the frequency of the transmitter:
Recommended distance			
Conducted interference induced by radio frequency electromagnetic fields according to IEC 61000-4-6 (GOST R 51317.4.6-99)	3B from 150kHz to 80MHz	V1 - 3 (V)	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$
Radio frequency electromagnetic field according to IEC 61000-4-3 (GOST 30804.4.3-2013)	3V/m from 80MHz to 2.5GHz	E1 - 3 (V / m)	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$ 80 MHz to 800 MHz
			$d = \left[\frac{7}{E_1} \right] \sqrt{P}$ 800 MHz to 2.5 GHz
Where P is the maximum transmitter output power in watts (W) according to the manufacturer's specifications, and the recommended distance in meters (m). d - recommended spatial separation, m; The field strength of radio waves from fixed radio transmitters, as determined by observations of the electromagnetic environment ^{a)} should be lower than the compliance level in each frequency band ^{b)} Interference may occur in the vicinity of equipment marked with :			

a) The field strength of radio waves from fixed radio transmitters, such as base stations of radiotelephone networks (cellular/wireless), and land mobile radio stations, amateur radio stations, AM and FM radio broadcast transmitters, television transmitters cannot be determined by calculation with sufficient accuracy. To do this, practical measurements of the field strength must be carried out. If the measured values at the location of the [ME EQUIPMENT or ME SYSTEM] exceed the applicable compliance levels, the operation of the [ME EQUIPMENT or ME SYSTEM] should be observed to verify that it is functioning properly. If monitoring reveals a deviation from normal operation, additional measures may need to be taken, such as reorienting or relocating the [ME EQUIPMENT or ME SYSTEM].

b) Outside the frequency range from 150 kHz to 80 MHz, the field strength should not exceed (V 1) V/m.

Table B.4 - Recommended values for spatial separation between portable and mobile radio frequency communications equipment and the system

Recommended separation distances between portable and mobile RF communications equipment and the system			
NOT INTENDED TO SUPPORT VITAL FUNCTIONS			
The system is intended for use in an electromagnetic environment in which radiated interference levels are controlled. The purchaser or user of the system can avoid electromagnetic interference by maintaining a minimum separation distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, taking into account the maximum output power of the communications equipment.			
Maximum transmitter output power (W)	Distance depending on transmitter frequency (m)		
	150 kHz ÷ 80 MHz	80 MHz ÷ 800 MHz	800 MHz ÷ 2.5 GHz
	$d = \left[\frac{3,5}{V_1} \right] \sqrt{P}$	$d = \left[\frac{3,5}{E_1} \right] \sqrt{P}$	$d = \left[\frac{7}{E_1} \right] \sqrt{P}$
0.01	0, 12	0, 12	0.23
0, 1	0.37	0.37	0.74
1	1, 17	1, 17	2, 33
10	3, 69	3, 69	7, 38
100	11, 67	11, 67	23, 33
Notes 1 At frequencies of 80 and 800 MHz, a higher field strength value is used. 2 The expressions given do not apply in all cases. The propagation of electromagnetic waves is affected by absorption or reflection from structures, objects and people. 3 When determining the recommended spacing values d for transmitters with a rated maximum output power not listed in the table, substitute the rated maximum output power P in watts specified in the transmitter manufacturer's documentation into the given expressions			

APPENDIX C

Weight and dimensions characteristics

Table C.1 – Weight and dimensional characteristics

Characteristic name	Meaning
Mouthpiece	
Length, mm, no more	111
Overall dimension in the transverse direction, mm, no more	16
Weight, g, no more	5
Mouthpiece according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes"	
Length, mm, no more	111
Overall dimension in the transverse direction, mm, no more	16
Weight, g, no more	5
Funnel mouthpiece	
Height, mm, no more	40
Overall dimension in the transverse direction, mm, no more	54
Weight, g, no more	5
Thermal paper	
Width, mm	56±2
External diameter, mm	23±2
Weight, g, no more	25
Battery pack	
LxWxH, mm, no more	65x44x20
Weight, g, no more	100
Analyzer power adapter	
LxWxH (body), mm, no more	95x55x90
Cable length, m, no more	3
Weight, g, no more	450
Analyzer on-board power adapter	
LxWxH (body), mm, no more	100x35x35
Cable length (for model with cable), m, no more	3
Weight, g, no more	100
Charger	
LxWxH, mm, no more	105x85x35
Weight, g, no more	100
AC adapter charger	
LxWxH (body), mm, no more	100x60x85
Cable length, m, no more	3
Weight, g, no more	300

End of Table C.1

On-board charger power adapter	
LxWxH (body), mm, no more	100x35x35
Cable length, m, no more	3
Weight, g, no more	150
USB cable	
Cable length, m, no more	3
Weight, g, no more	150
Analyzer power cable	
Cable length, m, no more	3
Weight, g, no more	150
Interface cable	
LxWxH (body), mm, no more	50x40x25
Cable length, m, no more	3
Weight, g, no more	150
Case (folded)	
LxWxH, mm, no more	180x80x50
Weight, g, no more	60
Case	
LxWxH, mm, no more	300x240x90
Weight, g, no more	600

APPENDIX D

Signs and symbols used in marking

Table D.1 – Signs and symbols used in marking

	Type approval mark
	Unified sign of product circulation on the market of the Eurasian Economic Union
	On/off button
	Product symbol type B
	C classification symbol according to the degree of protection against the penetration of water and solid particles
	D.C.
	Alternating current
	Non-ionizing radiation sign
	Manufacturer's name and address
	Prohibition on re-use
	Not sterile
	Do not use if packaging is damaged
	Carefully! Refer to instructions for use
	Refer to instructions for use
	Date of manufacture
	Batch number
	Don't burn

End of Table D.1

Mark of conformity for products distributed on the European market

Recycling symbol in accordance with Directive 2006/66/EC

Recycling symbol in accordance with Directive 2012/19/ EU

Recycle sign

Class II electrical product symbol

Product graphic symbol for indoor use

Positive polarity symbol

Information sign indicating the restriction of the content of hazardous substances

Manipulation sign "Fragile. Carefully". Fragility of the cargo. Careful handling of cargo

Manipulation sign "Keep away from moisture". The need to protect cargo from moisture

Manipulation sign "Top". Indicates the correct vertical position of the load

Manipulation sign "Keep away from fire." The need to protect cargo from fire.

Safety sign "Follow the operating instructions"

Serial number designation

Mini USB port designation

APPENDIX E

Maintenance Tools

Table E.1 – Maintenance means

Name and type of maintenance tool	Designation of a regulatory document regulating technical requirements and (or) metrological and basic technical characteristics of a maintenance device
Combination device * Testo 622	Temperature measurement range from + 10 °C to + 30 °C, limits of permissible absolute error ± 0.5 °C. Relative humidity measurement range from 30% to 80%, permissible absolute error limits $\pm 3\%$. Absolute pressure measurement range from 80 to 110 kPa, permissible absolute error limits ± 0.5 kPa.
Generator* of gas mixtures of ethanol vapor in air GUTH model 12 V 500 complete with standard samples of the composition of aqueous solutions of ethanol VRE-2 GSO 8789-2006	Limits of permissible absolute error in reproducing the mass concentration of ethanol in gas mixtures in the range (40 – 80) mg/m ³ ± 4 mg/m ³ . Limits of permissible relative error in reproducing the mass concentration of ethanol in gas mixtures in the range (over 80 – 2000) mg/m ³ $\pm 5\%$.
Rotameter* RM-1.6 GUZ	The upper limit of volume flow measurements is 1.6 m ³ /h; limits of permissible relative error $\pm 2.5\%$ of the upper measurement limit.
Digital multimeter * ARRA-72	The DC voltage measurement range is from 0 to 60 V, the DC current measurement range is from 0 to 6 A. The permissible absolute error limits for DC voltage measurements are $\pm (0.005 U_{\text{measured}} + 0.02)$ V. Limits of permissible absolute error in direct current measurements $\pm (0.01 \cdot I_{\text{measurement}} + 0.002)$ A.
Diaphragm compressor with flow regulator	Productivity not less than 10 l/min.
Distilled water	GOST R 58144-2018.

D.1 It is allowed to use other means not listed in Table D.1, but ensuring the required measurement accuracy.

D.2 All maintenance equipment marked * must be verified ; standard samples of the composition of aqueous ethanol solutions must have valid passports.

E.3 When carrying out periodic maintenance of analyzers, gas mixtures (hereinafter referred to as GS) are used with nominal values of the mass concentration of ethanol indicated in Table E.2.

Table E.2 – Metrological characteristics of HS

HS number (check points)	Nominal value of mass concentration of ethanol in HS, limits of permissible deviation, mg/l	Nominal value of the mass concentration of ethanol in aqueous solutions of ethanol ¹⁾ , limits of permissible deviation, mg/cm ³
1	0.00	distilled water
2	0.150 ± 0.015	0.386 ± 0.039
3	0.475 ± 0.048	1.22 ± 0.12
4	1.10 ± 0.11	2.83 ± 0.28

¹⁾ Standard samples of the composition of aqueous solutions of ethanol VRE-2 GSO 8789-2006. Relative error limits at P=0.95 : $\pm 1\%$.

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Technical Data Sheet

1 Basic information

Breath Alcohol Analyzer AJKOTEKTOP PRO-100 touch-M plus (hereinafter referred to as the analyzer) is designed for express measurements of the mass concentration of ethanol vapor in a selected sample of exhaled air.

contraindications for use or side effects.

Registration certificate No. _____ dated _____.

Number in the Federal Information Foundation for Ensuring the Uniformity of Measurements _____.

Manufacturer: Shenzhen Well Electric Co., Ltd., China.

Address: 1-3F, Building 101, No. 407, HeDongCun, HengKeng, GuanCheng Community, GuanHu Street, Longhua District, Shenzhen City, Guangdong, China, Tel: 86-755-83160728, Fax: 86-755-83160467, wellzp@well-co.com.

Authorized representative of the manufacturer (manufacturer): Limited Liability Company ALCOTEKTOR LLC (ALCOTEKTOR LLC), info@alcotector.ru.

Location address: 191036, St. Petersburg, 1st Sovetskaya st., 10, lit. Ah, placed. 2-N.

2 Metrological and technical characteristics

2.1 The metrological characteristics of the analyzer are given in Table 1.

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Measuring range mass concentration ethanol, mg/l	Limits of permissible error	
	absolute	relative
from 0 to 0.20 0 incl.	± 0.020 mg/l	-
St. 0.20 0 to 1.500	-	$\pm 10\%$

Note - The analyzer has a programmatically set minimum range of readings that are displayed on the analyzer screen and on paper in the form of zero readings: from 0.000 to 0.020 mg/l.

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2 If the analyzer readings exceed the upper reading limit, the message ">3,000 mg/l" is displayed on the analyzer screen and in the printed measurement protocol.

2.3 Price of the lowest digit of the scale, mg/l: 0.001.

2.4 There are no additional errors from the presence of unmeasured components.

2.5 Parameters of the analyzed gas mixture when a sample is supplied to the analyzer input (automatic sampling mode):

- flow rate of the analyzed gas mixture, l/min: at least 20.

- sample volume of the analyzed gas mixture, l: not less than 1.2.

2.6 Time to prepare for work after switching on, s:

- at ambient temperature St. +15 °C to +50 °C: no more than 5;

- at ambient temperature from -5 °C to +15 °C inclusive: no more than 10.

2.7 Measurement time after sampling a gas mixture with a mass concentration of ethanol of 0.000 mg/l (0.480 mg/l), s:

- at ambient temperature St. +15 °C to +50 °C:

- at ambient temperature from -5 °C to +15 °C inclusive:

2.8 Time to prepare for work after analyzing a gas mixture with a mass concentration of ethanol of 0.000 mg/l (0.480 mg/l), s:

- at ambient temperature St. +15 °C to +50 °C:

- at ambient temperature from -5 °C to +15 °C inclusive:

2.9 Operating time interval of the analyzer without adjusting the readings, months: no more than 5(10); no more than 5(30); not less than 12.

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2.10 The analyzer has built-in software TouchMplus.RU with the version number specified in 6.1.

2.11 A thermal paper roll is installed in the thermal printer, having the following overall dimensions: width 56±2 mm, outer diameter 23±2 mm.

2.12 Analyzer operating modes:

- measurement mode with data input;

- measurement mode without data input;

- screening mode;

- data transfer mode to PC;

- mode for adjusting and checking readings;

- touch screen calibration mode.

2.13 The analyzer has two modes of air sampling - automatic and manual.

2.14 The analyzer stores 10,000 results of the latest measurements in its memory.

2.15 Overall dimensions of the analyzer (LxWxH), mm, no more: 150x70x35.

2.16 Weight of the analyzer (with batteries, without thermal printer), g: no more than 280.

2.17 The analyzer can be powered electrically from:

- 4 AA batteries with a nominal voltage of 1.5 V (hereinafter referred to as batteries);

- 4 AA batteries with a voltage of 1.2 V and a capacity of at least 2200 mAh (hereinafter referred to as batteries);

- analyzer power adapter with an output voltage of ≈ 5.5 V and a current of at least 500 mA, connected to an AC mains;

- from the analyzer's on-board power adapter with an output voltage of ≈ 5 V and a current of at least 500 mA, connected to the vehicle's on-board network;

- from a USB port with a voltage of ≈ 5 V and a current of at least 500 mA via a USB cable, interface cable or analyzer power cable.

2.18 Number of measurements on the analyzer without replacing/charging batteries:

- without replacing batteries - at least 3000;

- without battery charge - at least 1500 (when using batteries with a capacity of at least 2200 mAh).

2.19 Terms of Use:

- ambient air temperature, °C: from -5 to +50;

- relative humidity

ambient air, %: from 10 to 100 (without moisture condensation);

- atmospheric pressure, kPa: from 84.0 to 106.7.

2.20 Service life of the electrochemical sensor installed in the analyzer, years: at least 2.

2.21 Average service life of the analyzer, years: 5.

2.22 Average time to failure, hours: 8000.

3 Individual characteristics

3.1 When sold, the analyzer has the following set of fields for data entry (in accordance with the analyzer operating manual):

☐ Option 1 ☐ Option 2 ☐ Other:

N o.	Name	Save*
1		
2		
3		
4		
5		
6		
7		
8		
9		
10		

*Data entered into the marked fields is stored in the analyzer's memory.

3.2 The analyzer has activated the function of adding signature fields to the measurement protocol:

☐ No ☐ Yes

3.3 The analyzer has the function of blocking the analyzer operation upon expiration of the verification period:

☐ No ☐ Yes

3.4 Upon request, the analyzer has the following functional features.

3.4.1 Provides the ability to work with the interface cable model:

 **It is not allowed to connect an interface cable to the analyzer whose model does not correspond to the one specified.**

3.4.2 Provides transmission of the measurement protocol:

☐ No ☐ Yes: ☐ to a thermal printer via Bluetooth
☐ to an external device via Bluetooth
☐ to an external device via an interface cable

3.4.3 Provides control of commands from an external device:

☐ No ☐ Yes: ☐ via Bluetooth
☐ via interface cable

3.4.4 Provides automatic switching on when power is supplied via the interface cable to the mini USB connector of the analyzer:

☐ No ☐ Yes

3.4.5 Provides data transfer to the PRO-100 touch-M plus Statistics program via a USB cable:

☐ No ☐ Yes

3.4.6 Provides determination of the coordinates of the measurement location:

☐ No ☐ Yes

4 Completeness

4.1 The analyzer delivery set is listed in Table 2.

Name	Designation	
Analyzer ethanol vapor in exhaled air	AJKOTEKTOP PRO-100 touch-M plus	
Mouthpiece	—	5 pieces
Funnel mouthpiece	—	1 PC. **
AA batteries 1.5 V	—	4 things.
USB cable	—	1 PC.
Case	—	1 PC.
Case	—	1 PC.
Thermal printer*	—	1 PC.
Battery pack*	—	1 PC. **
Thermal paper*	—	2 pcs. **
Charger AC adapter*	—	1 PC.
On-board charger power adapter*	—	1 PC.
Manual	RE 1.1	1 copy
Passport	—	1 copy

Notes:

- 1 The number of components in the table is indicated for the standard delivery set.
- 2 Components marked with "*" are supplied when ordering analyzers complete with a thermal printer.
- 3 The number of components marked with "***" may differ from those indicated in the table as agreed with the Customer.
- 4 The delivery set includes mouthpieces supplied by the manufacturer and/or mouthpieces in accordance with TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes".
- 5 Accessories supplied upon request and not included in the standard delivery package: charger, mains and on-board analyzer power adapters, AA 1.2 V batteries and battery charger, analyzer power cable, interface cable.
- 6 Components that are accessories: USB cable, funnel mouthpiece, mouthpieces according to TU 22.29.29-001-82139963-2017 (identical to TU 2291-001-82139963-2015), version "Alkotektor mouthpiece with two holes", thermal printer, battery pack, thermal paper, charger, mains and on-board power adapters for the analyzer, mains and on-board power adapters for the charger, power cable for the analyzer, interface cable.

5 Maintenance accounting

5.1 On-site maintenance includes:

- charging/replacing batteries at low voltage levels;
- charging the battery pack at low voltage (if applicable);
- refilling the thermal printer with thermal paper (if applicable);
- touch screen calibration - as needed;
- cleaning and disinfection - when soiled (frequency set by the user).

5.2 Periodic maintenance of the analyzer during the entire period of operation includes:

- technical condition monitoring - recommended before verification;
- checking the analyzer readings - it is recommended before verification (or more often in accordance with the frequency set by the user);
- adjusting the analyzer readings - in case of negative results of checking the readings or during verification as necessary.

5.3 It is recommended that records of maintenance and repairs performed, as well as records of changes made to the individual settings of the analyzer installed at the time

of sale according to 3.1, 3.2 and 3.3, and records of changes in the functional features of the analyzer according to 3.4, are noted in Table 3.

Note – Entries in the passport about work performed are made only by service center specialists and verifiers; If necessary, the user is recommended to keep a separate maintenance log for the analyzer.

6 Certificate of acceptance

6.1 Incoming control

Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M plus serial number _____ accepted in accordance with the current technical documentation and recognized as suitable for use.

Software version number _____

Date of manufacture _____

Representative of the Quality Control Department (authorized representative of the manufacturer):

initials, surname signature stamp Quality Control Department

6.2 Primary verification

Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M plus serial number _____

verified in accordance with document MP-242-2554-2023

“Breath Alcohol Analyzer AJIKOTEKTOP PRO-100 touch-M plus. Verification methodology”, agreed upon by the Federal State Unitary Enterprise “VNIIM named after. D.I. Mendelev” _____ 2023, in _____

organization, registration number of the accreditation certificate

and based on the results of the initial verification, it is recognized as complying with the metrological requirements established in the type description and suitable for use in the field of state regulation of ensuring the uniformity of measurements.

Verification date _____

Date of next verification no later than _____

Verifier: _____
signature initials, surname

Verification mark _____

The verification results are confirmed by information on the verification results entered into the Federal Information Fund for Ensuring the Uniformity of Measurements (FSIS “ARSHIN”).

Link
at FSIS
“ARSHIN”:

Т а б л и ц а 3 – Maintenance accounting

Initials, surname, signature									
Mark-stamp service center / stamp									
Completed work									
date									

Initials, surname, signature								
Service center stamp/stamp								
Completed work								
date								

7 Warranty obligations

7.1 ALKOTEKTOR LLC bears warranty obligations for the serial number _____ for 12 (twelve) months from the date of sale specified in 7.7.

7.7. When the analyzer is sold to the Buyer by the Seller, ALKOTEKTOR LLC is responsible for the quality and good condition of the analyzer for 12 (twelve) months from the date specified in 7.8, but not more than 15 (fifteen) months from the date specified in 7.7. In this case, the warranty obligations of ALKOTEKTOR LLC are valid only if there is the seal of ALKOTEKTOR LLC in 7.7 and the Seller's seal in 7.8. The implementation of warranty obligations is possible only upon presentation of this passport with the signature of the Buyer (operating organization) in 7.9.

Guaranteed shelf life is 6 months from the date specified in 7.7.

7.2 If any defects are found in the product during the warranty period, you must deliver the product to the ALKOTEKTOR LLC service center. If the claim is found to be justified, the Buyer (operating organization) or the Seller has the right to have the defects eliminated free of charge or replaced with a similar analyzer (at the discretion of ALKOTEKTOR LLC).

7.3 Acceptance of the product under warranty is carried out at the service center of ALKOTEKTOR LLC, tel. : (812) 320-22-97, at the address: 199178, St. Petersburg, emb. R. Smolenki, 5 -7, lit. A.

7.4 Procedure for accepting a product under warranty.

Fulfillment of warranty obligations is carried out in relation to analyzers that have a passport with a valid warranty, unbroken seals or brand marks (if any on the product) and a claim certificate. If the original passport is missing, ALKOTEKTOR LLC will remove the product from warranty service.

It is recommended to present the product in its entirety (excluding consumables).

The warranty obligations are fulfilled in relation to analyzers without mechanical damage, as well as other malfunctions resulting from:

- a) violation of the rules of operation of the product: use of the product for other purposes or not in accordance with the operating instructions for the analyzer,
- b) use of abnormal power supply,
- c) unqualified intervention, repairs carried out by unauthorized persons,
- d) transportation,
- e) accidents, natural disasters, fire, flooding, lightning, etc.

7.5 Transportation costs for delivery of a faulty product to the service center of ALKOTEKTOR LLC are at the expense of the Buyer (operating organization) or the Seller. The costs of sending the product to the Buyer after warranty repairs are at the expense of ALKOTEKTOR LLC.

7.6 If lost, the original passport cannot be restored.

7.7 date sales analyzer :

M.P.

_____ date of sale signature resp. face decoding signature

7.8 date implementation analyzer Seller :

M.P.

_____ Seller date signature resp. faces with transcript

7.9 I have read the warranty obligations:

_____ Buyer (operating organization) signature resp. person transcript signature date

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