Operating manual

Drying and heating ovens Binder in versions

Execution options:

E28 | **Natural convection ovens**

ED23

ED53

ED115

ED240

ED400

ED720

FD23 | **Forced convection ovens**

FD53

FD115

FD240

FED53 | **Natural and forced convection ovens**

FED115

FED240

FED400

FED720

Binder GmbH

- ▶ Address: Im Mittleren Oesch 5, 78532 Tuttlingen, Germany
- ▶ Post office box 102, 78502 Tuttlingen, Germany
- ▶ Tel.: +49 7462 2005 0 ▶ Fax: +49 7462 2005 100
- ▶ Internet: <http://www.binder-world.com> ▶ E-mail: info@binder-world.com
- ▶ Service - Hotline: +49 7462 2005 555
- ▶ Service - Fax: +49 7462 2005 93 555
- ▶ Service - E-mail: service@binder-world.com
- ▶ Service - Hotline, Russia and CIS: +7 495 988 15 16

Table of contents

DEAR CUSTOMERS!	5
1. SAFETY	11
1.1 Personnel qualification	11
1.2 Operating Manual	11
1.3 Legal considerations	12
1.3.1 Intellectual property	12
1.4 Structure of the safety instructions	12
1.4.1 Signal word panel	13
1.4.2 Safety alert symbol	13
1.4.3 Pictograms	13
1.4.4 Word message panel structure	14
1.5 Localization/position of safety labels on the device	14
1.6 Type plate	15
1.7 Transport marking	24
1.8 Basic instructions for the safe installation and use of ovens	25
1.9 Intended use	28
1.10 Foreseeable Misuse	29
1.11 Residual Risks	30
2. DESCRIPTION OF THE MEDICAL DEVICE	30
2.1 Medical device overview	33
2.1.1 E28 ovens	33
2.1.2 ED23, FD23, ED400, FED400 ovens	34
2.1.3 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens	34
2.1.4 Main power switch (ED720, FED720)	34
2.2 Triangular instrument panel	35
2.2.1 ED23, FD23, ED400, FED400 ovens	35
2.2.2 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens	35
2.3 Description of the oven options	36
2.3.1 Standard door or door with an inspection window	36
2.3.2 Door handle or door handle with a lock and a key	37
2.3.3 Chrome-plated racks	37
2.3.4 Stainless steel racks	38
2.3.5 Perforated stainless steel racks	38
2.3.6 Acoustic signaling	38
2.3.7 Independent thermal protection system, class 1	39
2.3.8 Independent thermal protection system, class 2 or class 3.1	39
2.3.9 Access port with silicone plug	39
2.3.10 Analog output 4-20 mA	39
2.3.11 HEPA filter (for FD53, FD115, FD240, FED53, FED115, FED240, FED400, FED720)	40
2.3.12 Inert gas connection system (for FD53, FD115, FD240, FED53, FED115, FED240 ovens)	40
2.3.13 Detachable mains cable	41
2.3.14 Temperature measurement calibration certificate (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary)	41
2.3.15 Calibration certificate for measuring spatial temperature (at 9 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)	49
2.3.16 Calibration certificate for measuring spatial temperature (at 27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)	56
3. COMPLETENESS OF DELIVERY, TRANSPORTATION, STORAGE AND INSTALLATION	64
3.1 Unpacking, checking equipment and completeness of delivery	64
3.2 Guidelines for safe lifting and transportation	64
3.3 Storage	65
3.4 Determining the installation location	65
3.5 Operating conditions	67
4. INSTALLATION	67

4.1	Mounting the tilt protection holders (chambers with "door with window" option)	67
4.2	Electrical connection	68
4.3	Ventilation pipe	69
4.4	Installing and loading racks	70
5.	COMMISSIONING	71
5.1	E28 ovens	71
5.1.1	Switching on the E28 oven with a timer 0-120 minutes	71
5.1.2	Setting the temperature	71
5.1.3	Fresh air supply	71
5.2	ED23, FD23, ED400, FED400 ovens	72
5.2.1	Turning on the oven	72
5.2.2	Heating function display	72
5.2.3	Air exchange	73
5.3	ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens	73
5.3.1	Turning on the oven	73
5.3.2	Functioning of the equipment when the door is opened	73
6.	BASIC SETTINGS OF THE CONTROLLER	73
6.1	ED23, FD23, ED400, FED400 ovens	73
6.1.1	Display / input of setpoints for temperature and ventilation (without setting the slow / ramp function)	73
6.1.2	Display / input of setpoints for temperature and ventilation (with setting of slow / ramp function)	74
6.1.3	Timing functions: Continuous operation and Operation under timer control	75
6.1.4	Custom settings	79
6.1.5	General recommendations	84
6.2	ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens	85
6.2.1	Overview and general settings on the R4 controller	85
6.2.2	Normal display mode	85
6.2.3	Setting the menu language	86
6.2.4	Setting date and time	87
6.2.5	Selecting the temperature unit	88
6.2.6	Set-point entry for temperature and fan speed	89
6.2.7	Adjusting the air flap position	91
6.2.8	Changing the passwords for user level and general controller functions	92
6.2.9	Slow/ramp function	94
6.2.10	Timer functions	98
7.	OVERTEMPERATURE PROTECTION (INDEPENDENT THERMAL PROTECTION SYSTEM)	109
7.1	Independent thermal protection system for E28 ovens	109
7.2	Independent thermal protection system for ED23, FD23, ED400, FED400 ovens	109
7.2.1	Independent thermal protection system Class 2	109
7.2.2	Independent thermal protection System class 3.1	111
7.3	Independent thermal protection system for ovens ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720	112
7.3.1	Independent thermal protection system Class 1	112
7.3.2	Independent thermal protection system class 2 ("temperature limiter")	113
7.3.3	Independent thermal protection system class 3.1 ("temperature protection" - protection against overtemperature)	113
7.3.4	Setting the safety controller set-point	114
7.3.5	Alarm message and proceeding in case of an alarm	114
8.	DATA RECORDING VIA USB INTERFACE FOR ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 OVENS	115
8.1	Starting data recording	115
8.2	Terminating data recording	117
9.	NETWORK CONFIGURATION FOR FED53, FED115, FED240, FED720 OVENS WITH ETHERNET INTERFACE	117
10.	CLEANING AND DECONTAMINATION	120
10.1	Cleaning oven	120
10.2	Decontamination	122

11. MAINTENANCE, TROUBLESHOOTING.....	123
11.1 Maintenance intervals	123
11.2 Troubleshooting	124
11.3 Sending the chamber back to the manufacturer / distributor / service center	125
12. END-OF-LIFE DISPOSAL.....	126
12.1 Disposal of the transport packing	126
12.2 Decommissioning	126
12.3 Disposal of the device on the territory of the Russian Federation	127
13. TECHNICAL DESCRIPTION	127
13.1 Factory calibration and adjustment	127
13.2 Definition of usable volume.....	127
13.3 Over current protection	128
13.4 Technical data for natural convection ovens	128
13.5 Technical characteristics of ovens with forced convection.....	131
13.6 Technical data for natural and forced convection ovens.....	132
13.7 Overview of oven equipment and functions.....	134
13.8 Accessories, spare parts and additional services (extract) - articles	135
13.8.1 Spare parts.....	135
13.8.2 Additional services from BINDER.	136
13.9 Medical device dimensions.....	137
13.9.1 Dimensions of oven E28 (mm).....	137
13.9.2 Dimensions of ovens ED23, FD23 (mm).....	138
13.9.3 Dimensions of ovens ED53, FD53, FED53 (mm).....	139
13.9.4 Dimensions of ovens ED115, FD115, FED115 (mm).....	140
13.9.5 Dimensions of ovens ED240, FD240, FED240 (mm).....	141
13.9.6 Dimensions of ovens ED400, FED400 (mm).....	142
13.9.7 Dimensions of ovens ED720, FED720 (mm).....	143
14. CERTIFICATES AND DECLARATIONS OF CONFORMITY.....	144
14.1 Registration certificate for medical device	150
15. PRODUCT REGISTRATION	152
16. ENVIRONMENTAL CLEANLINESS CERTIFICATE	163
17. MANUFACTURER'S WARRANTIES	155
18. INFORMATION ABOUT THE AUTHORIZED REPRESENTATIVE ON THE TERRITORY OF THE RUSSIAN FEDERATION.....	156

Dear customers!

The manufacturer, Binder GmbH company (hereinafter referred to as BINDER) presents the medical device: **Drying and heating ovens Binder in versions.**

I. Drying and heating ovens Binder with natural convection:

1. Model E28 including:

- 1.1 Drying and heating oven – 1 pc.
- 1.2 Rack, chrome-plated – from 2 to 4 pcs.
- 1.3 Rack, perforated, stainless steel – from 1 to 4 pcs. (if necessary)
- 1.4 Independent adjustable temperature safety device, class 1 – 1 pc. (if necessary)
- 1.5 Operating Manual – 1 pc.

2. Models ED23 including:

- 2.1 Drying and heating oven – 1 pc.
- 2.2 Standard door or Door with viewing window – 1 pc.
- 2.3 Door handle or Door handle with lock and key – 1 pc.
- 2.4 Chrome shelf – from 2 to 4 pcs.
- 2.5 Stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 2.6 Perforated stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 2.7 Acoustic alarm – 1 pc. (if necessary)
- 2.8 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 2.9 Access port with silicone plug – 1 pc. (if necessary)
- 2.10 Analog output 4-20 mA – 1 pc. (if necessary)
- 2.11 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 2.12 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 2.13 Operating manual – 1 pc.

3. Model ED53 including:

- 3.1 Drying and heating oven – 1 pc.
- 3.2 Standard door or Door with viewing window – 1 pc.
- 3.3 Door handle or Door handle with lock and key – 1 pc.
- 3.4 Chrome shelf – from 2 to 4 pcs.
- 3.5 Stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 3.6 Perforated stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 3.7 Acoustic alarm – 1 pc. (if necessary)
- 3.8 Independent thermal protection system class 1 – 1 pc.
- 3.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 3.10 Access port with silicone plug – 1 pc. (if necessary)
- 3.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 3.12 Removable network cable – 1 pc.
- 3.13 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 3.14 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 3.15 Operating manual – 1 pc.

4. Model ED115 including:

- 4.1 Drying and heating oven – 1 pc.
- 4.2 Standard door or Door with viewing window – 1 pc.
- 4.3 Door handle or Door handle with lock and key – 1 pc.
- 4.4 Chrome shelf – from 2 to 5 pcs.
- 4.5 Stainless steel shelf – from 1 to 5 pcs. (if necessary)
- 4.6 Perforated stainless steel shelf – from 1 to 5 pcs. (if necessary)

- 4.7 Acoustic alarm – 1 pc. (if necessary)
- 4.8 Independent thermal protection system class 1 – 1 pc.
- 4.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 4.10 Access port with silicone plug – 1 pc. (if necessary)
- 4.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 4.12 Removable network cable – 1 pc.
- 4.13 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 4.14 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 4.15 Operating manual – 1 pc.

5. Model ED240 including:

- 5.1 Drying and heating oven – 1 pc.
- 5.2 Standard door or Door with viewing window – 1 pc.
- 5.3 Door handle or Door handle with lock and key – 1 pc.
- 5.4 Chrome shelf – from 2 to 8 pcs.
- 5.5 Stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 5.6 Perforated stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 5.7 Acoustic alarm – 1 pc. (if necessary)
- 5.8 Independent thermal protection system class 1 – 1 pc.
- 5.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 5.10 Access port with silicone plug – 1 pc. (if necessary)
- 5.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 5.12 Removable network cable – 1 pc.
- 5.13 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 5.14 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 5.15 Operating manual – 1 pc.

6. Model ED400 including:

- 6.1 Drying and heating oven – 1 pc.
- 6.2 Standard door or Door with viewing window – 1 pc.
- 6.3 Door handle or Door handle with lock and key – 1 pc.
- 6.4 Chrome shelf – from 2 to 10 pcs.
- 6.5 Stainless steel shelf – from 1 to 10 pcs. (if necessary)
- 6.6 Perforated stainless steel shelf – from 1 to 10 pcs. (if necessary)
- 6.7 Acoustic alarm – 1 pc. (if necessary)
- 6.8 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 6.9 Access port with silicone plug – 1 pc. (if necessary)
- 6.10 Analog output 4-20 mA – 1 pc. (if necessary)
- 6.11 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 6.12 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 6.13 Operating manual – 1 pc.

7. Model ED720 including:

- 7.1 Drying and heating oven – 1 pc.
- 7.2 Standard door or Door with viewing window – 1 pc.
- 7.3 Door handle or Door handle with lock and key – 1 pc.
- 7.4 Chrome shelf – from 2 to 16 pcs.
- 7.5 Stainless steel shelf – from 1 to 16 pcs. (if necessary)
- 7.6 Perforated stainless steel shelf – from 1 to 16 pcs. (if necessary)
- 7.7 Acoustic alarm – 1 pc. (if necessary)
- 7.8 Independent thermal protection system class 1 – 1 pc.
- 7.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 7.10 Access port with silicone plug – 1 pc. (if necessary)

- 7.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 7.12 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 7.13 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 7.14 Operating manual – 1 pc.

II. Drying and heating ovens with forced convection

1. Model FD23 including:

- 1.1 Drying and heating oven – 1 pc.
- 1.2 Standard door or Door with viewing window – 1 pc.
- 1.3 Door handle or Door handle with lock and key – 1 pc.
- 1.4 Chrome shelf – from 2 to 4 pcs.
- 1.5 Stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 1.6 Perforated stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 1.7 Acoustic alarm – 1 pc. (if necessary)
- 1.8 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 1.9 Access port with silicone plug – 1 pc. (if necessary)
- 1.10 Analog output 4-20 mA – 1 pc. (if necessary)
- 1.11 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 1.12 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 1.13 Operating manual – 1 pc.

2. Model FD53 including:

- 2.1 Drying and heating oven – 1 pc.
- 2.2 Standard door or Door with viewing window – 1 pc.
- 2.3 Door handle or Door handle with lock and key – 1 pc.
- 2.4 Chrome shelf – from 2 to 4 pcs.
- 2.5 Stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 2.6 Perforated stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 2.7 Acoustic alarm – 1 pc. (if necessary)
- 2.8 Independent thermal protection system class 1 – 1 pc.
- 2.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 2.10 Access port with silicone plug – 1 pc. (if necessary)
- 2.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 2.12 HEPA filter – 1 pc. (if necessary)
- 2.13 System for connecting inert gas – 1 pc. (if necessary)
- 2.14 Removable network cable – 1 pc.
- 2.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 2.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 2.17 Operating manual – 1 pc.

3. Model FD115 including:

- 3.1 Drying and heating oven – 1 pc.
- 3.2 Standard door or Door with viewing window – 1 pc.
- 3.3 Door handle or Door handle with lock and key – 1 pc.
- 3.4 Chrome shelf – from 2 to 5 pcs.
- 3.5 Stainless steel shelf – from 1 to 5 pcs. (if necessary)
- 3.6 Perforated stainless steel shelf – from 1 to 5 pcs. (if necessary)
- 3.7 Acoustic alarm – 1 pc. (if necessary)
- 3.8 Independent thermal protection system class 1 – 1 pc.
- 3.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 3.10 Access port with silicone plug – 1 pc. (if necessary)
- 3.11 Analog output 4-20 mA – 1 pc. (if necessary)

- 3.12 HEPA filter – 1 pc. (if necessary)
- 3.13 System for connecting inert gas – 1 pc. (if necessary)
- 3.14 Removable network cable – 1 pc.
- 3.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 3.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 3.17 Operating manual – 1 pc.

4. Model FD240 including:

- 4.1 Drying and heating oven – 1 pc.
- 4.2 Standard door or Door with viewing window – 1 pc.
- 4.3 Door handle or Door handle with lock and key – 1 pc.
- 4.4 Chrome shelf – from 2 to 8 pcs.
- 4.5 Stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 4.6 Perforated stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 4.7 Acoustic alarm – 1 pc. (if necessary)
- 4.8 Independent thermal protection system class 1 – 1 pc.
- 4.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 4.10 Access port with silicone plug – 1 pc. (if necessary)
- 4.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 4.12 HEPA filter – 1 pc. (if necessary)
- 4.13 System for connecting inert gas – 1 pc. (if necessary)
- 4.14 Removable network cable – 1 pc.
- 4.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 4.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 4.17 Operation manual – 1 pc.

III. Drying and heating ovens with natural and forced convection

1. Model FED53 including:

- 1.1 Drying and heating oven – 1 pc.
- 1.2 Standard door or Door with viewing window – 1 pc.
- 1.3 Door handle or Door handle with lock and key – 1 pc.
- 1.4 Chrome shelf – from 2 to 4 pcs.
- 1.5 Stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 1.6 Perforated stainless steel shelf – from 1 to 4 pcs. (if necessary)
- 1.7 Acoustic alarm – 1 pc. (if necessary)
- 1.8 Independent thermal protection system class 1 – 1 pc.
- 1.9 Independent thermal protection system class 2 or 3.1 – 1 pc.
- 1.10 Access port with silicone plug – 1 pc. (if necessary)
- 1.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 1.12 Removable network cable – 1 pc.
- 1.13 HEPA filter – 1 pc. (if necessary)
- 1.14 System for connecting inert gas – 1 pc. (if necessary)
- 1.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 1.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 1.17 Operating manual – 1 pc.

2. Model FED115 including:

- 2.1 Drying and heating oven – 1 pc.
- 2.2 Standard door or Door with viewing window – 1 pc.
- 2.3 Door handle or Door handle with lock and key – 1 pc.
- 2.4 Chrome shelf – from 2 to 5 pcs.
- 2.5 Stainless steel shelf – from 1 to 5 pcs. (if necessary)

- 2.6 Perforated stainless steel shelf – from 1 to 5 pcs. (if necessary)
- 2.7 Acoustic alarm – 1 pc. (if necessary)
- 2.8 Independent thermal protection system class 1 – 1 pc.
- 2.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 2.10 Access port with silicone plug – 1 pc. (if necessary)
- 2.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 2.12 Removable network cable – 1 pc.
- 2.13 HEPA filter – 1 pc. (if necessary)
- 2.14 System for connecting inert gas – 1 pc. (if necessary)
- 2.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 2.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 2.17 Operating manual – 1 pc.

3. Model FED240 including:

- 3.1 Drying and heating oven – 1 pc.
- 3.2 Standard door or Door with viewing window – 1 pc.
- 3.3 Door handle or Door handle with lock and key – 1 pc.
- 3.4 Chrome shelf – from 2 to 8 pcs.
- 3.5 Stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 3.6 Perforated stainless steel shelf – from 1 to 8 pcs. (if necessary)
- 3.7 Acoustic alarm – 1 pc. (if necessary)
- 3.8 Independent thermal protection system class 1 – 1 pc.
- 3.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 3.10 Access port with silicone plug – 1 pc. (if necessary)
- 3.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 3.12 Removable network cable – 1 pc.
- 3.13 HEPA filter – 1 pc. (if necessary)
- 3.14 System for connecting inert gas – 1 pc. (if necessary)
- 3.15 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 3.16 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 3.17 Operating manual – 1 pc.

4. Model FED400 including:

- 4.1 Drying and heating oven – 1 pc.
- 4.2 Standard door or Door with viewing window – 1 pc.
- 4.3 Door handle or Door handle with lock and key – 1 pc.
- 4.4 Chrome shelf – from 2 to 10 pcs.
- 4.5 Stainless steel shelf – from 1 to 10 pcs. (if necessary)
- 4.6 Perforated stainless steel shelf – from 1 to 10 pcs. (if necessary)
- 4.7 Acoustic alarm – 1 pc. (if necessary)
- 4.8 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 4.9 Access port with silicone plug – 1 pc. (if necessary)
- 4.10 Analog output 4-20 mA – 1 pc. (if necessary)
- 4.11 HEPA filter – 1 pc. (if necessary)
- 4.12 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 4.13 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 4.14 Operating manual – 1 pc.

5. Model FED720 including:

- 5.1 Drying and heating oven – 1 pc.
- 5.2 Standard door or Door with viewing window – 1 pc.
- 5.3 Door handle or Door handle with lock and key – 1 pc.
- 5.4 Chrome shelf – from 2 to 16 pcs.

- 5.5 Stainless steel shelf – from 1 to 16 pcs. (if necessary)
- 5.6 Perforated stainless steel shelf – from 1 to 16 pcs. (if necessary)
- 5.7 Acoustic alarm – 1 pc. (if necessary)
- 5.8 Independent thermal protection system class 1 – 1 pc.
- 5.9 Independent thermal protection system class 2 or class 3.1 – 1 pc.
- 5.10 Access port with silicone plug – 1 pc. (if necessary)
- 5.11 Analog output 4-20 mA – 1 pc. (if necessary)
- 5.12 HEPA filter – 1 pc. (if necessary)
- 5.13 Calibration certificate for temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- 5.14 Calibration certificate for spatial temperature measurement (at 9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)
- 5.15 Operating manual – 1 pc.

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For proper use of the medical devices *Drying and heating ovens Binder in versions (hereinafter referred to as Ovens)*, it is necessary to read this Operating Manual and to observe its provisions carefully.

1. Safety

1.1 Personnel qualification

Personnel requirements

The oven should only be installed, set up, started, operated, cleaned and taken out of service by personnel who have the professional training, knowledge of relevant standards and experience to evaluate and perform their job responsibilities and to identify any potential hazards when their implementation. Such people must be sufficiently trained, instructed to perform work on the oven, and also have special permission to do this.

	Potential users The oven should only be operated by qualified laboratory or medical facility personnel who are specially trained for such purposes and who are familiar with the safety precautions.
	ATTENTION! • Installation, startup, operation, maintenance and repair of the oven by untrained, insufficiently qualified and unauthorized personnel is prohibited. • The human factor (lack of experience, qualifications, stress, fatigue, negligence) can cause an unfavorable event.

1.2 Operating Manual

Depending on the application and location of the cabinet, personnel must be provided with appropriate safety information when working with the cabinet in accordance with the Operating Instructions. The manual must be readable and written in the language of the employees working with the cabinet.

	 WARNING Failure to follow safety instructions. Serious damage and failure of the device. - Ø Read the instructions for safe handling of the device in this operating manual. - Ø Read the operating manual for the ovens carefully.
---	--

The instruction manual is part of the delivery set and should be kept in an easily accessible place in the immediate vicinity of the cabinet. Only specially trained personnel should be allowed to operate the device. It is important that all persons working with the cabinet read and understand the Operating Instructions. To avoid injury or damage, you must follow the safety instructions and the correct or correct procedure for operating the cabinet.

If you have any questions, we recommend contacting the manufacturer, BINDER, to clarify the data and applicability of these operating instructions.

	 DANGER
Danger arising from non-observance of instructions and safety requirements. Risk of serious injury to people or damage to the oven. Death risk. - Ø Follow the safety instructions in this manual. - Ø Follow the procedures in this manual. - Ø Before installing and operating the oven, you should carefully and completely read the manual for its correct operation. - Ø The operating manual must always be in a place accessible to personnel.	

	Ensure that all employees working with the oven read and fully understand this operating manual.
--	---

1.3 Legal considerations

This operating manual contains all the information necessary for the correct use, correct installation, start-up, operation and maintenance of the device.

Understanding and studying the provisions of this manual is a prerequisite for trouble-free operation and safety during operation and maintenance.

This operating manual cannot cover all conceivable applications. If you need more information or have specific problems that are not covered in this manual, please seek advice from your dealer or contact us directly.

Furthermore, we emphasize that the contents of this operating manual are not part of an earlier or existing agreement, description, or legal relationship, nor do they modify such a relationship. All obligations on the part of BINDER arise from the corresponding agreement, which also contains complete and correct information about the warranty. The statements in this manual neither augment nor restrict the contractual warranty provisions.

1.3.1 Intellectual property

Trademark information: All BINDER trade marks relating to products or services as well as the trade names, logos and product names used on the website, products and documents of BINDER company are trademarks or registered trademarks of BINDER (including BINDER GmbH, BINDER Inc.) in the U.S. and other countries and communities. Trademarks include word marks, position marks, word/figure marks, design configurations, and design patents.

Patent information: BINDER products, product categories and optional equipment may be covered by one or more patents and/or utility models in the U.S. and other countries and communities of states. Additional patent applications may also be pending in the U.S. and other countries and communities of states.

Please visit our website www.binder-world.com for more information.

1.4 Structure of the safety instructions

In this operating manual, the following safety definitions and symbols indicate dangerous situations following the harmonization of ISO 3864-2 and ANSI Z535.6.

1.4.1 Signal word panel

Depending on the probability and severity of consequences, potential dangers are identified with a signal word, the corresponding safety color and, if appropriate, the safety alert symbol.

 DANGER
Indicates an imminently hazardous situation that, if not avoided, will result in death or serious (irreversible) injury.

 WARNING
Indicates a potentially hazardous situation which, if not avoided, could result in death or serious (irreversible) injury

 CAUTION
Indicates a potentially hazardous situation which, if not avoided, may result in moderate or minor (reversible) injury

NOTICE
Indicates a potentially hazardous situation which, if not avoided, may result in damage to the product and/or impairment of its performance, as well as damage to nearby property.

1.4.2 Safety alert symbol

	Use of the safety alert symbol indicates a risk of injury. Observe all measures that are marked with the safety alert symbol in order to avoid death or injury.
---	---

1.4.3 Pictograms

Warning signs			
 Electrical hazard	 Hot surface	 Explosive atmosphere	 Stability hazard
 Lifting hazard	 Risk of corrosion and / or chemical burns	 Suffocation hazard	 Harmful substances
 Biohazard	 Pollution Hazard		

Mandatory action signs			
			
Mandatory regulation	Read operating manual	Disconnect the power plug	Lift with several persons
			
Lift with mechanical assistance	Environment protection	Wear protective gloves	Wear safety goggles
Prohibition signs			
			
Do NOT touch	Do NOT spray with water		
	Information to be observed in order to ensure optimum function of the product.		

1.4.4 Word message panel structure

Type / cause of hazard.
Possible consequences.
Instruction how to avoid the hazard: prohibition.
Instruction how to avoid the hazard: mandatory action.

Observe all other notes and information not necessarily emphasized in the same way, in order to avoid disruptions that could result in direct or indirect injury or property damage.

1.5 Localization/position of safety labels on the device

The following labels are located on the device:

Pictograms (Warning signs)	Service label
 Hot surface - E, ED, FD, FED: outer chamber door - ED, FD, FED: on chamber rear next to the exhaust duct - E: on the chamber top next to the air exhaust slots	
 Read operating manual	


Figure: Position of labels on the chamber front (example: ED, FD)

Keep safety labels complete and legible.


Replace safety labels that are no longer legible. Contact BINDER Service for these replacements.

1.6 Type plate

The type plate is located on the left-hand side of the chamber, bottom right-hand.

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Serial No. 00000000000000

Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
E 28		Model
Drying / heating oven		Oven name: Drying and heating oven
PY №ФС3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	00000000000000	Product Serial Number
Nominal temperature	230 °C 446 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	1.0	Overheating protection device class
Art. No.	9010-0001	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction

Indications of the plate (example)	Information
0,80 kW	Nominal power (kW)
3,5A	Nominal current
230 V / 50 Hz	Nominal voltage (V) $\pm$ 10% at the indicated power frequency (Hz)
230 V / 60 Hz	
1 N ~	Current type

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ED 23

Serial No. 00000000000000
 Made in Germany

Indications on the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
ED 23		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0358	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
0,80 kW		Nominal power (kW)
3,5A		Nominal current
230 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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
ED 400

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
ED 400		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0075	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
3,40 kW		Nominal power (kW)
7,4A		Nominal current
400 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
400 V / 60 Hz		
3 N ~		Current type

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
FD 23

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FD 23		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФЦ3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0360	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
0,80 kW		Nominal power (kW)
3,5A		Nominal current
230 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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FED 400

Serial No. 00000000000000
Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FED 400		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФЦ3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	00000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007

Indications of the plate (example)		Information
Class	2.0	Overheating protection device class
Art. No.	9010-0361	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
3,40 kW		Nominal power (kW)
7,4A		Nominal current
400 V / 50 Hz		Nominal voltage (V) $\pm 10\%$ at the indicated power frequency (Hz)
400 V / 60 Hz		
3 N ~		Current type

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Serial No. 00000000000000
Made in Germany

Drying and heating oven

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FD 53		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФС3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	00000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0346	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
1,10 kW		Nominal power (kW)
4,8A		Nominal current
230 V / 50 Hz		Nominal voltage (V) $\pm 10\%$ at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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
FD 115

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FD 115		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0354	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
1,30 kW		Nominal power (kW)
5,7A		Nominal current
230 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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

FD 240

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FD 240		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0347	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
2,30 kW		Nominal power (kW)
10,0A		Nominal current
230 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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Serial No. 00000000000000

Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FED 53		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	00000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature

Indications of the plate (example)		Information
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0348	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
1,10 kW		Nominal power (kW)
4,8A		Nominal current
230 V / 50 Hz		Nominal voltage (V) $\pm$ 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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Serial No. 00000000000000

Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FED 115		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФЦ3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	00000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0355	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
1,30 kW		Nominal power (kW)
5,7A		Nominal current
230 V / 50 Hz		Nominal voltage (V) $\pm$ 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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Serial No. 00000000000000

Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FED 240		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0349	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
2,30 kW		Nominal power (kW)
10,0A		Nominal current
230 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
230 V / 60 Hz		
1 N ~		Current type

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Serial No. 00000000000000

Made in Germany

Indications of the plate (example)		Information
BINDER		Manufacturer: BINDER GmbH Address: Im Mittleren Ösch 5, 78532 Tuttlingen / Germany Website: www.binder-world.com
FED 720		Model
Drying and heating oven		Oven name: Drying and heating oven
PY №ФC3 2010/08132 от 26.10.2010		Registration certificate for medical devices with number and date
Production date	xx-20xx	Date of manufacture of the product, month-year
Serial No.	000000000000000	Product Serial Number
Nominal temperature	300 °C 572 °F	Maximum temperature
IP protection	IP20	IP the type of protection in accordance with EN 60529
Safety device	DIN 12880	Overheating protection in accordance with DIN 12880:2007
Class	2.0	Overheating protection device class
Art. No.	9010-0356	Article number of the product
Project No.	---	Special application according to Project No.... (if applicable)
Built	20xx	Year of construction
4,50 kW		Nominal power (kW)
13,5A		Nominal current
400 V / 50 Hz		Nominal voltage (V) ± 10% at the indicated power frequency (Hz)
400 V / 60 Hz		
3 N ~		Current type

Symbols on the type plate	Information
	CE conformity marking
	The chamber is certified according to Customs Union Technical Regulation (CU TR) for the Eurasian Economic Union (Russia, Belarus, Armenia, Kazakhstan, Kyrgyzstan).
	Electrical and electronic equipment manufactured / placed on the market in the EU after 13 August 2005 is to be disposed of in a separate collection according to Directive 2012/19/EU on waste electrical and electronic equipment (WEEE).

1.7 Transport marking

Symbols used on the package labeling	Meaning
BINDER	Manufacturer
Made in Germany	Place of production
Model Art.No Serial No 1N ~	Product identification (label)
	Caution fragile goods
	Top of the product. Indicates the correct vertical position of the load

Symbols used on the package labeling	Meaning
	Keep away from moisture
	The maximum number of identical loads that can be stacked one on top of the other, where n is the limit
	Code for recycling corrugated cardboard for reuse
	Wood fiberboard recycling symbol
	Cardboard recycling and reuse symbol
	Symbol of comprehensive utilization and recycling of all transport and outer packaging of paper and cardboard. RECY is an association of manufacturers of corrugated cardboard and of paper recycling companies.

1.8 Basic instructions for the safe installation and use of ovens

Regarding the operation of the ovens and the determination of the installation location, please observe the safety rules and internal regulations in the laboratory room.

BINDER GmbH is only responsible for the safety of its equipment if experienced electricians or qualified personnel authorized by BINDER carry out all maintenance and repairs, and if all components related to the safety of the oven are replaced with original spare parts in the event of failure.

To operate the device, use only original BINDER accessories or accessories from third-party suppliers authorized by BINDER. The user is responsible for any risk caused by using unauthorized accessories.

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NOTICE	
	Danger of overheating due to lack of ventilation. Damage to the chamber. - Ø Do NOT install the chamber in unventilated recesses. - Ø Ensure ventilation for heat dissipation. - Ø Observe the prescribed minimum distances when installing the device (chap. 3.4)

Do NOT install or operate the ovens in hazardous locations.

DANGER	
	Danger of explosion due to combustible dusts or explosive mixtures in the vicinity of the oven. Serious injury or death from burns and / or explosion pressure. Ø Do NOT allow presence of any combustible dust or air and solvent mixtures near the oven.

The ovens do not dispose of any means of explosion protection.

	 DANGER
	Danger of explosion due to introduction of flammable or explosive substances in the device. Serious injury or death from burns and/or explosion pressure. Ø Do NOT introduce any substance into the oven which is combustible or explosive at working temperature in the inner chamber. Ø Do NOT introduce any combustible mixture in the inner chamber.

Any solvent contained in the charging material must not be explosive or inflammable. I.e., irrespective of the solvent concentration in the steam room, NO explosive mixture with air must form. The temperature inside the chamber must lie below the flash point or below the sublimation point of the charging material. Familiarize yourself with the physical and chemical properties of the charging material, as well as the contained moisture constituent and its behavior with the addition of heat energy.

Familiarize yourself with any potential health risks caused by the charging material, by the contained moisture constituent or by reaction products that may arise during the temperature process. Take adequate measures to exclude such risks prior to putting the chamber into operation.

	 DANGER
	Electrical hazard by water entering the chamber. Deadly electric shock. Ø The chamber must NOT become wet during operation, cleaning, or maintenance. Ø Do NOT install the oven in areas with high humidity and poor sealing. Set up the oven in a splash-proof manner.

The ovens are produced in accordance with VDE regulations and carefully tested in accordance to VDE 0411-1 (IEC 61010-1).

During operation and shortly after it, the temperature of the inner surfaces almost equals to the set-point. The inner chamber, exhaust duct, door window (in the presence of), and the door gaskets will become hot during operation of the oven, while the temperature of the outer surfaces of the oven does not exceed the permissible values.

Case temperature °C	Left-hand side	Front door	Right side	Upper side	Backside
E28	30	34	30	35	26
ED/FD 23	42	44	42	57	40
ED/FD/FED 53	34 -36	38	33-39	48-50	33-34
ED/FD/FED 115	32-33	43-48	31-33	45	31-34
ED/FD/FED 240	45	39-45	41-42	49	30
ED/FED 400	44-45	39-43	42-43	44	24-28
ED/FED 720	41-51	65-68	37-54,5	42-45,5	31-35

	CAUTION Interior chamber, ventilation duct, door window (optional) and door gaskets become hot during use. Burns. - Ø Do NOT touch the inner surfaces, exhaust duct, door window, access ports, door gaskets, or the charging material during operation. - Ø FD, FED: Do not place the power cable over the door gap when the chamber is hot after operation.
---	--

Requirements for the oven loading materials / substances/ products

None of the components of the loading material should be explosive, highly inflammable or a source of toxic gases. A mixture of any component of the load with air must NOT be explosive. The operating temperature must lie below the flash point or below the sublimation point of the loading material. The operating temperature must lie below the flash point or below the decomposition point of the loading material.

The loading material shall not contain any corrosive ingredients that may damage the machine components made of stainless steel, aluminum, and copper. Such ingredients include in particular acids and halides. Any corrosive damage caused by such ingredients is excluded from liability by BINDER GmbH.

It is NOT recommended to use the device for drying with the release of large amounts of steam, including leading to condensation.

	Based on the instructions in this operating manual, regular maintenance (chap. 11) is an integral part of correct operation.
	WARNING: If the customer has to use BINDER ovens outside the control of a service engineer, we strongly recommend storing non-recoverable material samples in at least two ovens, if possible.
	The charging material must not contain any corrosive components that can damage the oven parts. In particular, acids and halides. Any corrosion damage caused by such ingredients is excluded from the area of responsibility of BINDER GmbH.

The ovens do not dispose of any means of explosion protection.

	DANGER Explosion or implosion hazard. Danger of poisoning. Danger to life. - Ø Do NOT introduce any substance combustible or explosive at working temperature into the chamber, in particular no energy sources such as batteries or lithium-ion batteries. - Ø Explosive dust or air-solvent mixtures in the inner chamber are NOT ALLOWED. - Ø Do NOT introduce any substance which could lead to release of toxic gases.
---	---

Contamination of the ovens by toxic, infectious or radioactive substances must be prevented.

 	 WARNING Danger of intoxication and infection through contamination of the oven with toxic, infectious or radioactive substances. Damages to health. - Ø Protect the interior of the oven from contamination by toxic, infectious or radioactive substances. - Ø Take suitable protective measures when introducing and removing toxic, infectious or radioactive substances.
--	---

In case of foreseeable use of the device there is no risk for the user through the integration of the ovens into systems or by special environmental or operating conditions in the sense of EN 61010-1:2010. For this, the intended use of the oven and all its connections must be observed.

Connect only external devices to the oven Ethernet interfaces and USB which are compliant with the standards EN 61010-1:2010 or EN 60950-1:2006 modified.

Products made of heat-resistant materials can be placed in the oven, for example, surgical, dental, gynecological instruments, glass and glass instruments, heat-resistant syringes and needles, as well as other medical products that can be sterilized with hot air both in the package and without it.

1.9 Intended use

Purpose

The ovens are designed for drying, heat treatment and sterilization of solid or pulverized substances and materials as well as products made of heat-resistant materials using heat supply.

Are used

The ovens are applicable in organizations carrying out medical activities or laboratory practice.

Indications for use

Not applicable for this product.

Contraindications and side effects

Not applicable for this product.

Oven parameters/operating modes

Ovens can be used, for example, for drying glassware, as well as for warm storage of liquids in containers. They are not suitable for drying processes where large amounts of steam are generated leading to condensation. It is necessary to comply with the requirements for the installation location (chapter 3.4) and the requirements for loading (chapter 1.8) of the ovens described in the Operation Manual. Only trained personnel (Chapter 1.1) who are familiar with the Operation Manual are allowed to work with ovens.

Ovens can provide sterilization modes at temperatures of 85/120 and 180 °C.

Sterilization at 180 degrees (recommended): Set the set point on the controller to 190 degrees to ensure that each internal surface reaches 180 degrees. It is necessary to allow the camera to warm up and then transfer it to stable mode for an hour. Then let it run for at least 30 minutes for effective sterilization. Then turn off the oven and wait until it cools down. This temperature also ensures that bacterial spores are inactivated.

Sterilization at 85 degrees: Set the set point on the controller to 95 degrees and allow the oven to heat up. When the actual value becomes equal to the set value, let the oven operate for at least 30 minutes. Then turn it off and wait until it cools down.

Sterilization at 120 degrees: Set the set point on the controller to 130 degrees and allow the oven to heat up. When the actual value becomes equal to the set value, operate the oven for at least 30 minutes. Then turn it off and wait until it cools down.

The consumer can also use ovens for heat treatment and drying.

Depending on the oven model, the technical specifications indicate the maximum temperature settings. For all ovens the maximum temperature is 300°C, with the exception of E28, where it is 230°C.

E28 ovens: select the temperature mode and set the mechanical timer for the selected duration, maximum - 120 minutes, or set it to infinity (symbol ∞), so that the chamber works continuously until the timer is manually turned off.

ED23, FD23, ED400, FED400 ovens: set the desired temperature and select on the controller between "Continuous operation" or "Operation by timer". Maximum timer operating time: 99 hours 59 minutes. In "Continuous Operation" mode, the cabinet operates continuously until it is manually turned off.

ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens: set the desired temperature and select on the controller between "Continuous operation" or "Timer operation". Maximum timer operating time: 9 days 23 hours 59 minutes. Without the timer function, the cabinet runs continuously until it is manually turned off.

1.10 Foreseeable Misuse

Unacceptable actions that entail consequences and pose a risk to the equipment, regardless of its inherently safe design:

- Failure to comply with the operating instructions
- Failure to comply with information and warnings on the cabinet (eg control device messages, safety labels and warning signals)
- Failure or delay in performing maintenance and testing.
- Ignores signs of wear and damage
- Use of non-recommended materials.
- Failure to comply with technological parameters for processing certain materials
- Installation, testing, maintenance or repair in the presence of solvents
- Installation of spare parts or use of additional devices and technologies that are not approved by the manufacturer.
- Non-use or modification of protective systems, operation of the cabinet without original protective devices
- Ignoring messages about cleaning and disinfecting the cabinet
- Spillage of water or cleaning agent into the camera, or water entering the camera during operation, cleaning or maintenance.
- Perform cleaning operations while the cabinet is running
- Operating a cabinet with a damaged housing or power cable
- Continued operation of the cabinet if a malfunction is detected
- Closing the air outlet in the upper part of the housing (E28) and cabinet dampers with foreign objects, especially metal ones

To prevent the described and other risks of incorrect operation, it is recommended that appropriate work instructions and standard operating procedures (SOPs) be prepared for users.

1.11 Residual Risks

Oven risk assessments and analysis were performed which are recorded and maintained in the Risk Analysis file. The risk analysis addressed potential hazards that are common when the medical device is used as intended, as well as any hazards that can arise unexpectedly (for example, an explosion, a fall due to acceleration or deceleration, etc.), the risk was assessed and controls were implemented for its mitigation.

Raw materials for devices are purchased from approved suppliers. The device is manufactured in a controlled and verified shop environment following the approved manufacturing procedures and appropriate controls involving quality assurance.

The unavoidable design features of an oven/a chamber, as well as its proper field of application, can also pose risks, even during correct operation. These residual risks include hazards which, despite the inherently safe design, existing technical protective equipment, safety precautions and supplementary protective measures, cannot be ruled out.

Messages on the chamber and in the Operating Manual warn of residual risks. The consequences of these residual risks and the measures required to prevent them are listed in the Operating Manual. Moreover, the operator must take measures to minimize hazards from unavoidable residual risks. This includes, in particular, issuing operating instructions.

2. Description of the medical device

Drying and heating ovens (ovens) described in this Operating manual are manufactured by BINDER, a leading global supplier of specialized chambers for simulating environmental conditions in research and industrial laboratories.

The advantages of using BINDER ovens are: uniform temperature distribution; fast, uniform drying; wide temperature range; low power consumption; high performance and environmental safety. Ovens are designed for continuous operation and can operate continuously, with breaks only for maintenance.

The ovens are intended for drying, heat treatment and sterilization of solid or pulverized substances and materials as well as products made of heat-resistant materials using heat supply.

The E and ED natural convection ovens solve the daily tasks of drying and sterilizing at temperatures up to 230 °C (E28) and 300 °C (ED) and storing at a precisely controlled temperature. Thanks to natural convection of air with a high rate of air exchange, thermal processes take place with significantly higher efficiency.

FD forced convection ovens are used when fast drying and sterilization is required. Absolutely uniform temperature distribution, fast dynamic response and a specially designed air turbine allow significant time savings.

FED natural and forced convection ovens simultaneously satisfy several consumer needs, including sterilization, drying and controlled storage in a wide range of temperature steps with high precision.

Ovens are heated electrically (heating elements are used). The air flow in the E28 and ED ovens occurs naturally by gravity (natural convection). Forced circulation is organized in FD and FED ovens by means of a fan.

The air change function ensures high accuracy in maintaining temperature and temperature transitions due to direct and uniform air circulation inside. In FD and FED ovens, fans ensure that the desired temperature is reached and maintained with the required accuracy, and also increase the rate of temperature recovery. An air flap in the exhaust duct (ED, FD, FED) or beneath the exhaust slots (E28) serves to adjust the air change.

The air flap in E28, ED23, FD23, ED400, FED400 ovens is opened/closed with a mechanical arm located on the oven controller or inside the oven (for E28 only).

The position of the air flap in ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens to adjust the fresh air entry through the exhaust duct is adjusted with buttons on the controller.

If necessary, the introduced fresh air is cleaned by means of a high efficiency HEPA class H14 filter acc. to EN 1822:2009. This fine filter is located on the left side of the oven and can be easily replaced. All ovens, with the exception of the E28, have a rear exhaust duct. E28 has exhaust holes at the top.

Ovens, depending on the model in accordance with DIN 12880:2007, are equipped with an independent thermal protection system of class 1 (overheating device), an independent thermal protection system of class 2 (temperature limiter), and also, upon customer request, an independent thermal protection system of class 3.1 (temperature protection). The purpose of the independent thermal protection system is to protect the ovens, its surroundings and the loaded material from unacceptable excessive temperatures (see Chapter 7).

Sterilization: The recommended sterilization temperature is 180°C. Sterilization time: for at least 30 minutes, starting from the moment the sterilization temperature is reached. Ovens can provide additional sterilization modes at temperatures of 85 and 120 °C. The consumer may also be guided by other sterilization rules, if available.

The ovens consist of a supporting frame, inside of which there is a working volume (rectangular chamber) with a layer of thermal insulation. Access to the working volume is through a door with a polyamide handle; the loading opening is located horizontally. Depending on the principle of loading and unloading the chamber, the oven is one-sided. At the customer's request, a standard door is installed in the oven at the factory, or a door with a viewing window can be installed instead of a standard door. Also, upon customer request, instead of a standard handle, a door handle with a lock and key can be installed at the factory. A door seal is fixed along the perimeter of the doorway. The door seal is made of heat-resistant silicone rubber. A hermetically sealed seal is guaranteed for the entire period of operation of the product.

Two or several shelves are installed inside the working volume. Thoughtful fastening of shelves to the side walls allows you to place products of sufficiently large mass on them without breaking or deforming the shelves themselves, as well as rearranging them to the required height level. When loading the oven, the correct orientation of the shelves, the maximum permissible load per shelf and the maximum permissible total load must be observed (see technical specifications).

Oven racks: The oven is equipped with chrome-plated racks as standard. Chrome-plated racks can be optionally replaced with stainless steel racks. Besides, ovens can be optionally equipped with perforated stainless steel racks. The number of racks can be optionally increased (see the technical data sections).

The inner chamber, the preheating chamber and the inner side of the doors are made of V2A stainless steel (material No. 1.4301 in Germany, USA equivalent to AISI 304). When the ED, FD and FED ovens are operated at temperatures above 150°C, exposure to oxygen in the air can cause discoloration of metal surfaces (brownish yellow or light blue) as a result of natural oxidation. Such staining is harmless and in no way affects the quality of the equipment functioning. The ovens body, all corners and edges are powder coated in RAL 7035. Shelves and other surfaces are made of stainless steel or chrome, which generally eliminates the risk of corrosion. All parts of the ovens are easy to clean, which helps protect the processed materials from contamination.

The ovens are manufactured according to VDE requirements and have been thoroughly tested according to VDE 0411-1 (IEC 61010-1). During operation and for a short time after its termination, the

temperature of the internal surfaces is almost equal to the set value. During operation of the oven, the internal chamber, ventilation pipe, door viewing window (if equipped), and door gaskets become hot, while the temperature of the external surfaces of the oven does not exceed permissible values.

Case temperature °C	Left-hand side	Front door	Right side	Upper side	Backside
E28	30	34	30	35	26
ED/FD 23	42	44	42	57	40
ED/FD/FED 53	34 -36	38	33-39	48-50	33-34
ED/FD/FED 115	32-33	43-48	31-33	45	31-34
ED/FD/FED 240	45	39-45	41-42	49	30
ED/FED 400	44-45	39-43	42-43	44	24-28
ED/FED 720	41-51	65-68	37-54,5	42-45,5	31-5

Depending on the temperature, the heating mode of ED oven “adjusts” to the conditions of door opening. In FD and FED ovens, when the door is opened, the heating and fan are turned off and do not turn on while the door remains open. The force when opening/closing the door and door handle does not exceed 150H.

Ovens E28 has a mechanical thermostat and timer.

The microprocessor controller will allow you to pre-select, install and program the drying/sterilization processes of the loaded material, and carry out natural or forced convection. All programming functions are clear and clearly organized. Temperature ranges are given in the technical specification’s sections.

FED400 ovens are equipped with RS 422 interface as a standard, and ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens have USB connectors; optionally FED53, FED115, FED240, FED720 can be equipped with Ethernet interface for computer connection.

ED720 and FED720 ovens are equipped with four castors. Both front castors can be locked by brakes.

At the consumer's request (if necessary), the standard configuration of the oven can be upgraded for ease of use, depending on the consumer's goals, with the following components:

- Stainless steel shelf
- Perforated stainless steel shelf
- Acoustic alarm
- Access port with a silicone plug
- Analog output 4-20 mA
- HEPA fresh-air filter (for FD53, FD115, FD240, FED53, FED115, FED240, FED400, FED720 ovens)
- System for inert gas connection (for FD53, FD115, FD240, FED53, FED115, FED240 ovens)

At the consumer's request (if necessary), in order to confirm the quality of the oven, the manufacturer can provide verification certificates complete with the ovens:


- Certificate of calibration of temperature measurement (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).
- Calibration certificate for spatial temperature measurement (9/27 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)

2.1 Medical device overview


2.1.1 E28 ovens

	(1) Control thermometer – allows you to control the temperature up to 250°C
	(2) Mechanical Timer 0-120 minutes and Main Switch at the same time
	In position zero (0), heating is turned off. In position (∞) heating is constantly on In position (10-120) heating is turned on until the set time
(3) Green indicator lamp (ready for operation) (4) Thermostat (for setting the temperature up to 230°C)	
	(5) Yellow indicator lamp (heating in progress)
	(6) Locking device for thermostat (8) Fuse (9) Power cable


2.1.2 ED23, FD23, ED400, FED400 ovens

2.1.3 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens


2.1.4 Main power switch (ED720, FED720)

2.2 Triangular instrument panel


2.2.1 ED23, FD23, ED400, FED400 ovens

- (1) Controller display
- (2) Controller function buttons
- (3) Shutdown indicator
- (4) Control knob for the system of independent thermal protection against overheating
- (5) Lever for ventilation flap



2.2.2 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens

- (7) Controller display
- (8) Controller function buttons
- (9) USB connector


2.3 Description of the oven options


2.3.1 Standard door or door with an inspection window

For ED23, FD23, ED400, FED400 ovens	For ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens
Standard door	
	
Door with an inspection window	
An inspection window in the door with an insulating glass unit is provided for easy control of the sample - observation of samples, process and sterilization indicators.	
	


2.3.2 Door handle or door handle with a lock and a key

For E28, ED23, FD23, ED400, FED400 ovens	For ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens
Door handle	
	
For ED23, FD23, ED400, FED400 ovens	For ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens
Door handle with a lock and a key	
A lock with a key can be installed on the oven door handle to limit the access of unauthorized persons to the oven chamber.	
	


2.3.3 Chrome-plated racks

	
General view	View in the chamber

2.3.4 Stainless steel racks

2.3.5 Perforated stainless steel racks

2.3.6 Acoustic signaling

The function of acoustic signaling of excess temperature can be installed in any oven, except for E28. This function is activated and deactivated on the dashboard in manual mode.

	Disabling the acoustic signaling does not affect the operation of the independent thermal protection system.
--	--

For ED23, FD23, ED400, FED400:

The buzzer can be turned on/off with a switch.

Position 0 = buzzer off (not activated)

Position 1 = buzzer on (activated)

If the buzzer is activated, when the temperature limit set in the independent thermal protection system of class 2 or class 3.1 is exceeded, in addition to the flashing of the red warning lamp on the dashboard, an acoustic signal is activated.

For ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720:

Position Off = buzzer off (not activated)

Position On = buzzer on (activated)

The buzzer can be turned on/off by pressing the OK button.

If the buzzer is activated, when the temperature limit set in the independent thermal protection system of class 2 or class 3.1 is exceeded, in addition to the flashing of the alarm message on the controller display, an acoustic signal is generated.

The buzzer can be turned on/off by pressing the OK button.


2.3.7 Independent thermal protection system, class 1

See Chapter 7.3.1

2.3.8 Independent thermal protection system, class 2 or class 3.1

See Chapter 7.2.1, 7.2.2, 7.3.2, 7.3.3

2.3.9 Access port with silicone plug

The access port with a silicone plug is a hole in the wall of the oven with a lined pipe insert for the input of external measuring devices, sensors in order to independently control the functional parameters of the oven, for example, temperature. It comes with a silicone plug.

2.3.10 Analog output 4-20 mA

Pin Arrangement in DIN Socket for analog output


The 4-20 mA analogue output option, installed in the oven, allows the transfer of temperature signal data to external data logging systems or devices.

The connection is made via a 6-pin DIN socket on the back of the chamber. DIN plug included.

The connection is made via the DIN socket on the back of the chamber. Suitable DIN connector included

Analog output 4-20 mA, DC	Pin 1: temperature - Pin 2: temperature +
Temperature range:	ED, FD, FED: 0°C to +300°C

2.3.11 HEPA filter (for FD53, FD115, FD240, FED53, FED115, FED240, FED400, FED720)

HEPA filter	Standard mounting location
	

With this option, the supply air is purified by means of a high-efficiency filter from particles of less than submicron size – HEPA class H14 (in accordance with EN 1822:2009). If necessary, you can replace the filter element by removing the metal filter cover on the left side of the device.



2.3.12 Inert gas connection system (for FD53, FD115, FD240, FED53, FED115, FED240 ovens)

The oven is equipped with two ports for inert gas (nitrogen or inert gases).

The ports are located in the center on the top of the unit and on the bottom right of the rear of the oven. Each of these ports can be used as an inlet or outlet, depending on the nature of the inert gas:

- lighter gas (nitrogen, helium): the top port serves as an inlet
- heavy gas (e.g. argon): the bottom port serves as an inlet

This distinction is important when the fan speed is reduced.

	
---	--

With this option, the oven is additionally sealed, resulting in less leakage. To avoid overpressure, the device is not completely sealed. The sealing reduces the amount of vapors passing through the oven body, which can be released when the feed is heated. Discharge of vapors through a ventilation duct, for example to a contaminated air collection device, further reduces unwanted emissions. In order to avoid overpressure, the device is not completely sealed.

The device is not completely sealed. Gases from the inside of the oven can be released into the surrounding atmosphere.


It is necessary to comply with the restrictions regarding the maximum permissible concentrations of emitted substances in the work area, as well as comply with all applicable standards.

Any harmful gases emitted must be removed by good ventilation in the work area or by using a suitable flue gas system. If necessary, install the device in the immediate vicinity of ventilation.

The air flap does not completely cover the ventilation duct of the device. The supplied plug prevents the release of vapors or loss of inert gas supplied through the vent. Due to special thermal resistance requirements, use only the original supplied plugs.

When drying, always remove the cap to prevent the scattering of escaping vapors, which can lead to condensation of vapors in the inner chamber of the device.

2.3.13 Detachable mains cable

The ED53, ED115, ED240, FD53, FD115, FD240, FED53, FED115, FED240, ovens are equipped with a detachable mains cable.

The mains cables of E28, ED23, ED400, ED720, FD23, FED400, FED720 ovens are not detachable and form part of the oven.

2.3.14 Temperature measurement calibration certificate (at 1 point) in accordance with the BINDER factory standard – 1 pc. (if necessary).

This is a paper document for additional confirmation of the temperature characteristics of the oven at a given temperature. One measurement point (usually the center of the oven) is checked for one temperature value.

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Measured values within the specification

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We hereby confirm that above mentioned chamber was calibrated following the quality assurance system according to DIN EN ISO 9001 (Certificate Reg. No. 12 100 27534 TMS)

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All technical data (Product catalog) are specified for units in standard model.

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This calibration certificate may not be reproduced other than in full except with the permission of the Binder GmbH and the Person in charge/

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This test certificate is only valid for the above-mentioned chamber with the controller and sensors that were used in this measurement.

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The measurements were performed after the unit is one hour in a thermally equilibrated condition.

The unit was measuring in empty condition with three racks.

Measuring points: in accordance to BINDER factory standard following DIN 12880.

Duration of measurement: At least one hour.

This test certificate is only valid for the named device under simultaneous use of controller and sensor.

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1 temperature sensor is placed in the unit in accordance to BINDER factory standard.

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*Limits of the technical data

All indicated data for series devices are typical average values. The maximum permissible deviation shall not exceed or fall below 20% of the specified average value.
This criterion concerns the specified technical data of the chambers and the creation of calibration certificates.

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**** Value not specified by BINDER**

Reason: No information given in data sheet, or modified chamber construction

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Applicable Documents:

You can order the BINDER Factory Standard. Please contact BINDER Service.

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2.3.15 Calibration certificate for measuring spatial temperature (at 9 points) in accordance with the BINDER factory standard – 1 pc. (if necessary)

This is a paper document for additional confirmation of the temperature characteristics of the oven at a given temperature. In addition to the center of the oven, additional points are measured at 1 temperature value.

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The measurements were performed after the unit is one hour in a thermally equilibrated condition. The unit was measuring in empty condition with three racks.

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3. Completeness of delivery, transportation, storage and installation

3.1 Unpacking, checking equipment and completeness of delivery

After unpacking, please check the chamber and its optional accessories, if any, based on the delivery receipt for completeness and for transportation damage. Inform the carrier immediately if transportation damage has occurred.

The final tests of the manufacturer may have caused traces of the racks on the inner surfaces. This has no impact on the function and performance of the chamber.

Please remove any transportation protection devices and adhesives in/on the chamber and on the doors and take out the operating manuals and accessory equipment.

When transporting in cold conditions, move the device to a warmer installation location; the camera will need time to warm up. Condensation may form. Before unpacking, you should wait at least one hour for the device to reach ambient temperature and ensure that the camera is completely dry before starting.

	 CAUTION Sliding or tilting of the chamber. Damage to the chamber. Risk of injuries by lifting heavy loads. - Ø Do NOT lift or transport the chambers/ovens using the door handle or the door. - Ø Do NOT lift ED720 and FED720 ovens by hand. ➤ Lift the 23,28, 53 and 115 size chambers from the pallet at four lower corners with the aid of 2 people, 240 size chambers with the aid of 4 people. ➤ Lift 720 size chambers from the pallet using technical devices (fork lifter). Set the fork lifter only from the rear in the middle of the chamber. Make sure to place all the lateral supports of the chamber on the forks.
--	---

If you need to return a chamber/an oven, please use the original packing and observe the guidelines for safe lifting and transportation (chap. 3.2).

For disposal of the transport packing, see chap. 12.1.

Note on second-hand chambers (Ex-Demo-Units):

Second-hand chambers are chambers that have been used for a short time for tests or exhibitions. They are thoroughly tested before resale. BINDER ensures that the chamber is technically sound and will work flawlessly.

Second-hand chambers are marked with a sticker on the chamber door. Please remove the sticker before commissioning the chamber.

3.2 Guidelines for safe lifting and transportation

The front castors of ED720, FED720 ovens can be blocked by brakes. Please move empty ovens with castors only when empty and on an even surface, otherwise the castors may be damaged. After operation, please observe the guidelines for temporarily decommissioning (chap. 12.2).

	 CAUTION Sliding or tilting of the chamber. Damage to the chamber. Risk of injuries by lifting heavy loads. ➤ Transport the chamber only in its original packaging. ➤ Secure the chamber with transport straps for transport. ⊗ Do NOT lift or transport the chamber using the door handle or the door. ⊗ Do NOT lift 720 size chambers by hand. ➤ Lift 23, 28, 53 and 115 size chambers from the pallet at four lower corners with the aid of 2 people, 240 size chamber with the aid of 4 people, and place it on a transport pallet with wheels. Push the pallet to the desired site and then lift the chamber from the pallet at its four lower corners. ➤ Place ED720, FED720 chamber using technical devices (fork lifter) on the transport pallet. Set the fork lifter only from the rear in the middle of the chamber. Make sure to place all the lateral supports of the chamber on the forks. ➤ Transport ED720, FED720 chambers ONLY with the original transport pallet. Set the fork lifter only to the pallet. Without the pallet the chamber is in imminent danger of overturning!
---	---

- Permissible ambient temperature range during transport: -50 °C to +60 °C.
- Permissible ambient humidity: maximum 70%, non-condensing.

You can order transport packing and pallets for transportation purposes from BINDER Service.

3.3 Storage

Intermediate storage of the chamber is possible in a closed and dry room. Observe the guidelines for temporary decommissioning (chap. 12.2).

- Permissible ambient temperature range during storage: -10 °C to +60 °C.
- Permissible ambient humidity: max. 70 % r.H., non-condensing

When after storage in a cold location you transfer the chamber to its warmer installation site, condensation may form. Before start-up, wait at least one hour until the chamber has attained ambient temperature and is completely dry.

3.4 Determining the installation location

The oven should be placed on a flat, non-vibrating surface in a well-ventilated, dry place and leveled using a level. The installation location must be able to support the weight of the device (see technical data, chapter 13.4-13.6). The cabinets are intended for indoor use only.

When installing and operating the cabinet, it is necessary to provide access to normal ventilation (window) or establish a ventilation system. When working with inert gases, the necessary placement condition includes installing an additional ventilation system or connecting the outlet to the ventilation pipe.

	NOTICE	
	Danger of overheating due to lack of ventilation. Damage to the chamber. - Ø Do NOT install the chamber in unventilated recesses. ➤ Ensure ventilation for heat dissipation. ➤ Observe the prescribed minimum distances when installing the device.	

Installation height: maximum 2000 m above sea level.

When installing several ovens of the same size in a row, it is necessary to leave a distance between them of at least 250 mm. The distance to the ceiling is at least 100 mm. Also, when installing, it is necessary to maintain the distance to the wall at the back and the distance to the wall at the side.

Model		E28	ED23 FD23	ED53 FD53 FED53	ED115 FD115 FED115	ED240 FD240 FED240	ED400 FED400	ED720 FED720
Distance to wall, rear (minimum)	mm	100	100	160	160	160	100	160
Distance to wall, side (minimum)	mm	160	100	100	100	100	160	100
Free space above the chamber	mm	200	--	--	--	100	100	100

Two devices with a volume of up to 115 l can be stacked on top of each other. To do this, place a soft rubber pad under all four feet of the upper unit to prevent slipping.

	NOTICE	
	Risk of equipment slipping. Risk of equipment damage. - Ø When stacking, place rubber pads under all four feet of the upper device. - Ø Stack only ovens of the same size.	

Ovens sizes 28, 240, 400 and 720 must NOT be stacked.

	NOTICE	
	Danger when stacking ovens Risk of equipment damage. Ø Do not stack the ovens	

To completely separate the device from the power supply, you must disconnect the power plug. Install the device in a way that the power plug is easily accessible and can be easily pulled in case of danger. Do not conduct the power cable above the exhaust duct.

There is no risk to the consumer in the event of short-term overvoltages covered by the standard EN 61010-1: 2010.

DO NOT install or operate the chamber in potentially explosive areas.

	 DANGER
	Danger of explosion.
	Danger to life.
	Ø Do NOT operate the chamber in potentially explosive areas. Ø KEEP explosive dust or air-solvent mixtures AWAY from the vicinity of the chamber

3.5 Operating conditions

Permissible ambient temperature range during operation: from +10°C to +40°C. At elevated ambient temperatures, temperature fluctuations may occur. Permissible ambient humidity: maximum 70%, non-condensing

	The ambient temperature must not be significantly higher than +25°C, with which the given technical data are related. In case of different environmental conditions, deviations from the given data are possible.
---	---

4. Installation

4.1 Mounting the tilt protection holders (ovens with door with window)

For chambers equipped with the “door with window” option it is recommended to install the supplied tilt protection.

Preparing the tilt protection holders

- Depending on the desired wall distance, you can bend the tilt protection holders accordingly.


Figure: Variable length of the tilt protection holder depending on the bend

Installation on the chamber (A)

- Plug the two tilt protection holders each with the tab on the provided spot on the edge of the rear panel. The screw holes in the rear wall and in the tilt protection holders must align.
- Fix the tilt protection holders each with one of the supplied screws on the chamber rear wall.

Wall mounting (B)

- Then fix both tilt protection holders on the wall, each with 2 screws Ø 6mm

4.2 Electrical connection

The ovens are supplied ready for connection. They come complete with a main connecting cable 1800/2000 mm long (+/- 50 mm) depending on the model. All ovens are equipped with an IEC plug.

Model	Mains connector	Rated voltage +/- 10% at a given frequency	Current type	Fuse
E 28	Grounding plug	230 V at 60 Hz	1N~	6.3 A
ED23 FD23	Grounding plug	230 V at 50 Hz 230 V at 60 Hz	1N~	10 A
ED400 FED 400	Grounding plug	400 V at 50 Hz 400 V at 60 Hz	3N~	---
ED 53 FD 53 FED 53	Grounding plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6.3 A
ED 115 FD 115 FED 115	Grounding plug	230 V at 50 Hz 230 V at 60 Hz	1N~	6.3 A
ED 240 FD 240 FED 240	Grounding plug	230 V at 50 Hz 230 V at 60 Hz	1N~	12.5 A
ED720 FED720	Grounding plug	400 V at 50 Hz 400 V at 60 Hz	3N~	---

- The outlet must also provide protective earthing. Make sure that the protective earthing of the socket and plug are compatible!

	 DANGER
	Danger of electric shock due to unconnected protective conductor. Lethal electric shock. ➤ Make sure that the camera plug and socket are matched and the electrical protection wires of the camera and home installation are securely connected.

- Use only original BINDER connecting cables.

	NOTICE
	Danger of incorrect power supply voltage. Damage to the equipment. Ø Check the power supply voltage before connection and start-up. Ø Compare the power supply voltage with the data indicated on the type plate.

- When connecting, you must comply with the rules set by your local energy supply company. We recommend using an RCD (residual current device).
- FD, FED: Do not place the power cable over the door gap when the chamber is hot after operation.

- Pollution degree (according to IEC 61010-1): **2**
- Overvoltage category (according to IEC 61010-1): **II**
- Electromagnetic compatibility (according to EN 61326-1.): class A

See also electrical data (Chapter 13.4-13.5).

	To completely separate the device from the power supply, you must disconnect the power plug. Install the device in a way that the power plug is easily accessible and can be easily pulled in case of danger.
---	---

4.3 Ventilation pipe

Active suction from the chamber must only be affected together with external air. Therefore, the chamber's exhaust air duct shall not be immediately connected to an active exhaust system.

When connecting to an active exhaust system, proceed as follows:

- Perforate the connecting piece between the exhaust air duct and the exhaust system.

Or

- Use an exhaust air funnel placed in a distance of 3-5 cm from the exhaust air duct. The funnel's opening must be at least twice as large as the diameter of the exhaust air duct.

	If improperly connected to an active exhaust/ventilation system, the spatial temperature exactitude (uniformity), the heating-up and recovering times as well as the maximum temperature of the chamber may be negatively affected.
--	---


The exhaust duct on the chamber rear (ED/FD/FED) and the exhaust slots on top of the housing (E 28) will become hot during operation.







	 CAUTION
	Danger of burning by touching hot chamber parts during operation. Burns. Ø ED/FD/FED: Do NOT touch the exhaust duct during operation. Ø E28: Do NOT touch the air exhaust slots on top of the housing during operation and do not cover them.

4.4 Installing and loading racks

Observe the correct orientation of the racks:

The lateral brackets must be above the rack surface when inserting the rack.

		
Standard chrome rack	Stainless steel rack	Perforated stainless steel rack
		
Chamber view of chrome and stainless steel racks	Chamber view of perforated stainless steel rack	

When loading the oven, observe the maximum permissible load per rack and the maximum permissible total load (see technical data, chap. 16).

Observe the correct orientation of the racks:

 	 WARNING
	Overloading of racks. Risk of injury. ⊘ Do not exceed the maximum permissible rack load. ⊘ Do not exceed the maximum permissible total load. ⊘ Use only racks designed for this oven model. ⊘ The correct orientation of the racks must be observed when installing them. ⊘ Place the load carefully on the racks.

5. Commissioning

Ovens may release odors in the first few days after commissioning. This is not a quality defect. To reduce odors quickly we recommend heating up the oven to its maximal temperature for one day and in a well-ventilated location.

5.1 E28 ovens

5.1.1 Switching on the E28 oven with a timer 0-120 minutes


Figure: Timer 0-120 minutes (E 28)

In position zero (0), heating is turned off. The wardrobe doesn't work.

When setting the time in the range from 10-120 (the timer knob must be turned clockwise), the heating is turned on - the cabinet operates until the set time. After the time has elapsed, heating will stop automatically. The green indicator lamp will go out.

In the position ∞ (infinity symbol) – continuous mode heating is switched on (heating occurs continuously). It is installed by turning it counterclockwise, while the green indicator light is on. The cabinet operates for a long time in continuous operation mode at the set temperature until the user turns it off.

5.1.2 Setting the temperature

Set the thermostat knob (4) to the desired temperature. Lock it by turning the thermostat lock (6).

The yellow control lamp (5) will light up, signaling the heating process. When the operating temperature is reached, the yellow indicator light will start flashing, signaling that the thermostat is operating.

To ensure accurate temperature control, always set the temperature by turning the thermostat knob (4) clockwise. Before setting any new temperature, turn the thermostat knob counterclockwise to the initial position (zero).

5.1.3 Fresh air supply

The ventilation flap inside the chamber on the ceiling serves to regulate the supply of fresh air.

If the ventilation flap is fully open, the spatial temperature accuracy is reduced to 5%.

The openings for exhausting air from the inner chamber become hot during the operation of the device, do not cover them.

	NOTICE
	Danger of overheating. Damage to the device. Ø DO NOT touch the air vents of the inner chamber while the oven is in operation. Ø DO NOT cover the inner chamber air vents located on the top cover of the oven.

5.2 ED23, FD23, ED400, FED400 ovens

5.2.1 Turning on the oven

Insert the plug into a suitable socket (chap. 4.2).

ED400, FED400: switch on the oven at the main power switch (10) (chap. 2.1.4).

The green "Standby" indicator is on


Keep  pressed until the display lights up.

The controller is now in normal display mode (display of current value).

If the oven is in operation (temporary functions "Continuous operation" or "Timer operation" with a set time that counts down, Chapter 6.1.3), the **actual temperature value** is displayed (example: 22°C).

If the oven controller is controlled by the time function "Timer operation" (no time programmed or the countdown has already expired (Chapter 6.1.3), the device is inactive (does not heat up). The display alternately shows the actual **temperature value** (for example: 22°C) and "tOff":

An independent thermal protection system must be installed in accordance with the selected temperature setting (chap. 7.2.1.).

5.2.2 Heating function display

Heating and fan (for FD23, FED400) are running when the red heating control indicator in the lower right corner of the display starts flashing slowly depending on the heat demand (example: 70°C)

5.2.3 Air exchange

Opening the air flap in the exhaust duct serves to adjust the air change.

Without connecting the suction pump to the ventilation pipe:

- In the ED series, the circulation of fresh air can be increased by means of a ventilation pipe. The ventilation flap in the exhaust duct is designed to regulate the flow of incoming air.
- In the FD and FED series, if the air flap is open and the fan is running, fresh air is drawn in through the ventilation gaps.
- If the air flap is completely open, the spatial temperature accuracy can be negatively influenced.


Figure: Air flap adjustment

5.3 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens

5.3.1 Turning on the oven

Insert the connector into the appropriate socket (chap. 4.2).

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If the controller display does not show anything other than the standby symbol, press and hold the standby button until the display lights up.

The controller display switches to normal mode (chap. 6.2.2). If the timer function was active before the oven was turned off, it will be shown on the controller display.

5.3.2 Functioning of the equipment when the door is opened

ED: Depending on the temperature, heating performance may be adapted when opening the door.

FD, FED: When the door is opened, the heating and fan are turned off and will not turn on as long as the door remains open.

6. Basic settings of the controller

6.1 ED23, FD23, ED400, FED400 ovens

6.1.1 Display / input of setpoints for temperature and ventilation (without setting the slow / ramp function)

The controller setting is identical on all ED23 FD23, ED400 devices. The temperature range of the controller is up to 300°C, the display accuracy is one degree.

The device is running; the controller is now in normal display mode (displaying the current value). Current temperature value (for example: 22°C) is displayed:

1. Click the  button

The display shows alternately "SP" and the previous **temperature setpoint** (example: 60°C):

2. Enter the value using the   buttons.

The desired setpoint can be selected within a temperature range of 5°C above room temperature and up to 300°C

3. Wait 2 seconds for the entered temperature value to be accepted by the controller (display will flash once).
4. For FED400: Press the  button to proceed to entering the fan speed.
5. For FED400: The display shows alternately "n" and the previous **setpoint for the fan speed** (example: 100%):

6. For FED400: Set the desired fan speed using the   buttons.

The fan speed can be set with a value from 0% - 100%.

7. For FED400: Wait 2 seconds for the entered value to be accepted by the controller (the display will flash once).
8. Press the  button to return to the normal data display (display of the current value) (or this will happen automatically after about 30 seconds).

If the value is changed, it is necessary to re-install the independent thermal protection system.

6.1.2 Display / input of setpoints for temperature and ventilation (with setting of slow / ramp function)

If the slow / ramp function was previously set (chap. 6.1.4.2):

Press the  button in normal display / temperature display mode (for FED400) during slow/ramp to display the change in current value according to the selected gradient in addition to the entered temperature value and fan speed (for FED400).

The device is running; the controller is now in normal display mode (displaying the current value). **The current temperature** is displayed (for example: 22°C):

1. Click the  button.

The display shows alternately "SPr" and the **current setpoint temperature**, which changes according to the selected gradient (for example: 42°C):


This ramp setpoint is only displayed, not adjusted.

- Click the  button.


The display shows alternately “SP” and the previous **temperature setpoint** (example: 60°C):


- Enter a setpoint between 0 and 300 using the   keys.


- Wait 2 seconds for the entered temperature value to be accepted by the controller (display will flash once).
- For FED400: Press the  button to proceed to entering the fan speed.
- For FED400: The display shows alternately “n” and the previous **setpoint for the fan speed** (example: 100%):

- For FED400: Set the desired fan speed using the   buttons.

- For FED400: Wait 2 seconds for the entered value to be accepted by the controller (the display will flash once).
- Press the  button to return to the normal data display (or this will happen automatically after about 30 seconds).

6.1.3 Timing functions: Continuous operation and Operation under timer control

Press the time control button . The timer shows its current timing function. There are two possible timing functions:

Continuous operation mode

The display shows alternately “t1” and the time function “Continuous operation” “t inf”

Heating occurs continuously, regardless of the timer setting.

Operation under timer control

The display shows alternately “t1” and the **time** until the end of the countdown or “tOFF”

 or		 Remaining time (example: 28 min) - The timer is counting down Heating and ventilation (for FD) continues until the time set on the timer expires (for FED400 - depending on the selected timer function).
		 The timer is not programmed or the countdown has ended "t off" If the timer expires, heating and ventilation (for FD23) stops (for FED400 - the behavior of the device depends on the previously selected timer function).

Press the  button to return to the normal data display (display of the current value) (or this will happen automatically after about 30 seconds).

6.1.3.1 Switching between Continuous operation and Operation under timer control

Press the time control button .

The controller timer shows its current time function. When the "Continuous operation" function is active, "t1" and "t inf" are displayed alternately. When the "Operation under timer control" function is active, "t1" is displayed alternately with the time remaining until the end of the operation or "tOff".


If the Timer is counting down in the timing function ("t1" is displayed alternately with the counting down time), the timer must first be set to zero ("0") (Chapter 6.1.3.3). Now "t1" is displayed alternately with "tOff" and the controller can be switched to the timing function "Continuous operation".

 Continuous operation		
 2 seconds  2 seconds		
 Operation under timer control		 The timer is not programmed or the countdown has ended


Press the  button to return to the normal data display (display of the current value) (or this will happen automatically after about 30 seconds).

6.1.3.2 Continuous operation mode

1. Press the time control button . The timer shows the current timing function.
2. If necessary, switch to the "Continuous operation" function with the  button.
The display shows alternately "t1" and the time function "Continuous operation"
"t inf" (for FED400 only) or "tOFF" for other models:

or

3. Press the  button to return to the normal data display (display of the current value) (or this will happen automatically after about 30 seconds).

Current temperature value (for example: 22°C) is displayed:


The controller now operates with the entered setpoints (Chapter 6.1.1) in "Continuous operation" mode. Heating and ventilation (for FD23) occurs continuously, regardless of the timer setting.

To cancel "Continuous operation", proceed accordingly:

1. Press the time control button .
2. Switch to "Operation under timer control" by pressing the  button for 2 seconds (chap. 6.1.3.3).

6.1.3.3 Operation under timer control

1. Press the time control button . The timer shows the current timing function.
2. If necessary, switch to "Operation under timer control" with the  button. The display shows alternately "t1" and the **time until the end of the countdown** or "tOFF"


3. Set the desired time [hh, mm] with the arrow keys  .

Maximum setting: 99 hours 59 minutes. (for FD23, ED23, ED400)

4. Wait 2 seconds for the entered temperature value to be accepted by the controller (display will flash once).

The display shows alternately “t1” and the **time** until the end of the countdown.

The time will start counting immediately after the entered value takes effect. Heating and ventilation (for FD) are active until the time set on the timer has elapsed: for the FED400, the activity of the device is determined by the previously selected timer function.

5. Press the  button to return to the normal data display (display of the current value) (or this will happen automatically after about 30 seconds).

The **current temperature value** is displayed (for example: 22°C):

The controller operates with the entered setpoints (Chapter 6.1.1) until the set time has expired. Heating and ventilation (for FD) are active until the time set on the timer has elapsed. Heating activity (for FED) depends on the time entered and the timer function selected in the user menu.

To check or change the remaining timer time, press the time control button  in the normal display (current value display).

The display shows alternately “t1” and the **time until the end of the countdown**

After the set time has elapsed, the display shows the current temperature value (example: 22°C) and “tOff”:

Heating is not active now. Ventilation continues to work (FED), or not active (FD)

6.1.3.4 Timer functions for FED400

The device provides three different timer functions:

- Switch-off delay (setting “0”):

After the set time has elapsed, the heating is turned off.

- Switch-off delay depending on temperature (setting “1”)

The set time starts to count only when the current temperature value reaches the set one minus 1°C. After the set time has elapsed, the heating is turned off.


- Switch-on delay (setting “2”)

After the set time has elapsed, the heating is switched on and the oven operates in “Continuous operation” mode.


1. Press and hold the  button for about 5 seconds.
The display alternately shows "unit" and the current temperature unit setting.

2. Press the  button again.
The display alternately shows "rASd" and the currently set gradient value:

3. Press the  button again.
The display alternately shows "Adr" and the currently set value of the device address:

4. Press the  button again.
The display alternately shows "tFCt" and the current set value of the **Timer functions**:

5. Set the desired value for the timer functions 0, 1 or 2 using the   buttons.
The set value will automatically take effect after 2 seconds.

6.1.4 Custom settings

By pressing the  button in the normal display mode (displaying the current value) and holding it for 5 seconds, you enter the user menu.

The settings in this menu affect the behavior of the controller.

User level overview:



Figure: Press and hold the  button for about 5 seconds.

Press the  button to return to the normal temperature setting display.

OR:

After approx. 30 seconds, the controller will automatically return to normal/live data display.

All installations can be carried out independently (as described in separate sections) or sequentially one after the other in one simple process.

The set parameters are not deleted when the main unit is turned off or in the event of a power failure.

6.1.4.1 Switching the temperature scale between Celsius °C and Fahrenheit °F

If necessary, the temperature display can be changed as follows:

1. Press and hold the  button for about 5 seconds.

The display alternates between “unit” and the current temperature unit setting.

2. Use the   buttons to set the desired unit.
3. The set unit is automatically adopted after 2 seconds.

	C = degrees Celsius 0°C = 31°F Recalculation: F = degrees Fahrenheit 100°C = 212°F [Value in °F] = [Value in °C] * 1.8 + 32
---	---

When entering values for the slow / ramp function (Chapter 6.1.4.2), the set unit will be used as the basis.

	If the unit is changed, the set temperatures and limits are converted accordingly.
---	--

6.1.4.2 Slow/ramp function

You can program temperature ramps in order to extend heating up times. This may be necessary in some cases to prevent thermal stress in the material during heating. Temperature ramps should only be used as needed. Using temperature ranges may result in slower heating times. When the slow/ramp function is turned off, the chamber will heat up with its maximum heating capacity.

Entering °C/min or °F/min denotes the nominal gradient value and limits the maximum temperature rise to this value. Through heating and evaporation, the energy acquired by the material can lead to lower temperature gradients.

Slow/ramp heating is performed from the previously entered value to the new setpoint. The temperature must be adjusted to the starting setpoint.

Enter the settings in 3 steps:

1. Enter the set value for the start temperature of slow/ramp heating. Allow the temperature to adjust to this value.
2. Set the function to the desired gradient.

The gradient can be set between 1.0°C/min and 10.0°C/min.

If the gradient is set to 0, then the slow/ramp function is disabled, i.e. the heating rate will be as high as possible.

If the gradient is specified by a different number, for example, 3, then the heating rate will be as close as possible to 3°C/min.

A heating rate of 4°C/min is considered the most realistic maximum possible.

3. Enter the setpoint (final temperature).

The range should only be set if necessary. A setting of “0” means that the slow/ramp function is disabled. The device will heat up with the highest possible heat input.

1. Press and hold the  button for about 5 seconds.

The display alternates between “unit” and the current **temperature unit setting**.

2. Press the  button again.

The display alternately shows “rASd” and the **currently set gradient value**:

3. Set the desired range gradient using the   buttons (setting the gradient in °F or °C according to the setting in the Chapter 6.1.4.1).
4. The set value is automatically activated after 2 seconds.
During slow/ramp heating, the set value (SPr) is continuously increased according to the entered gradient from the previously entered to the newly entered (SP).

For a description of the display when setting the slow/ramp function, see Chapter 6.1.2.

6.1.4.3 Oven addressing

If several ovens (optional) are networked and connected to a PC using the software, each device must be assigned a unique address.

Addresses are set on the controller as follows:

1. Press and hold the  button for about 5 seconds.

The display alternates between “unit” and the current **temperature unit setting**.

2. Press the  button again.

The display alternately shows “rASd” and the **currently set gradient value**:

3. Press the  button again.

The display alternately shows “Adr” and the currently set value of the **device address**.

4. Set the desired value with the   buttons.

You can enter an address value between 1 and 30.

4. The set value will automatically take effect after 2 seconds.

6.1.4.4 Setting the interface mode and printer interval

1. Press the  button for about 5 seconds

The display alternates between “unit” and the current **temperature unit setting**.

2. Press the  button again.

The display alternately shows “rASd” and the **currently set gradient value**:

3. Press the  button again.

The display alternately shows “Adr” and the currently set value of the **device address**.

4. Press the  button again.

The display alternately shows “tFCt” and the current set value of the **timer function**:


5. Press the  button again.

The display will alternately show “PFct” and the current **interface mode setting**:

6. Set the desired interface mode with the   buttons.

Installations: Modbus = “0” printer = “1”

The set value will automatically take effect after 2 seconds.

If interface mode “1” (printer) has been selected, the printer interval for automatic output can be set in an additional menu step.

7. Press the  button again.

The display alternately shows “Prt” and at the input level the current setting of the **printer interval**:

8. Set the desired value with a number from the range 0 - 255 using the   buttons.

The printer intervals via the RS 422 interface can be set between 1 and 255. Setting “0” indicates that printer intervals are off.

The **printer** records temperature data at a set interval.

The set value will automatically take effect after 2 seconds

6.1.4.5 Temperature programming example

The device should warm up to a temperature of 50°C, maintain this temperature for 3 hours and then turn off.


1. During normal display status, hold down the  button for about 5 seconds and then press a few more times until “tFCt” is displayed.
 - Select timer function “1” = “Temperature dependent switch-off delay function” (chapter 6.1.4.4)
2. Press the  button during normal display status.
 - Enter the fixed value “50” (Chapter 5.2.1.).
3. During normal display status, press the time control button . The controller shows the current time function
 - If necessary, select the time function “Timer operation” (Chapter 6.1.3.1)
 - In the entry level, enter the desired time “3.00” (chap. 6.1.3.)

6.1.5 General recommendations

- | | |
|---|--|
|  | After approx. 30 seconds after the last entry, the controller returns to the normal display state (displaying the current value). |
|  | The entry of the function setpoints (Chapter 6.1.1), temporary functions (Chapter 6.1.3) and the call up of the user menu (Chapter 6.1.3.4) can only be selected from the normal display state. |
|  | During the entry of setpoints, time functions and selection of functions in the user menu, the corresponding  or  button must be pressed within 1 second. Pressing the button for a shorter time interval will be ignored by the controller. |
|  | After a power outage (for example, in case of a power drop), the timer will return to its previous status. The remaining time, if any, will continue to count down. |
|  | If the value is changed, it is necessary to re-install the independent thermal protection system. |

6.2 ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens

6.2.1 Overview and general settings on the R4 controller

Buttons on the controller

		The arrow buttons serve to navigate and to enter the values
		The OK button serves to select the parameters and to confirm the entered values
		The Back button serves to reach the preceding level
		If the Standby button is pressed down for approx. 3 seconds, the display changes to standby mode. To activate the display, press down the standby button again for approx. 3 seconds
		Display in standby mode with standby symbol

Status symbols on the controller display


	Heating active		Timer operation		Overtemperature alarm of the safety controller		Standby mode
---	----------------	---	-----------------	---	--	---	--------------

6.2.2 Normal display mode



	Normal display with chambers without fan (ED) or with fixed fan speed (FD)
---	--

	Normal display with chambers with adjustable fan speed (FED)
---	--

6.2.3 Setting the menu language

	From Normal display
without fan 5x with fan 6x	 with the arrow-up button to the user menu
	Confirm with OK.
	Enter the password (factory setting: 00 00) and confirm each entry with OK.
4 x	 with the arrow-up button to the language setting menu.
	The current menu language is shown.
	Press OK to select the menu language.
	The setting flashes.
	Select the setting with the arrow buttons
	and confirm with OK.
2x	 Back to Normal display.

There are the following options:

German: 	English: 	French: 
--	---	---

6.2.4 Setting date and time


	From Normal display
without fan 5x with fan 6x	 with the arrow-up button to the user menu
	Confirm with OK.
	Enter the password (factory setting: 00 00) and confirm each entry with OK.
	The current date is shown.

	Press OK to set the year .
	The setting flashes.
	Enter the year with the arrow buttons (any setting)
	and confirm with OK.

	Press OK to set the month .
	The setting flashes.
	Enter the month with arrow buttons (1 to 12)
	and confirm with OK.

	Without the optional real time clock, these settings must be repeated when the power supply is interrupted.
---	---

	Press OK to set the day .
---	----------------------------------

	The setting flashes.
	Enter the day with arrow buttons (1 to 31)
	and confirm with OK.

	Press OK to set the hour .
	The setting flashes.
	Enter the hour with arrow buttons (0 to 23)
	and confirm with OK.

	Press OK to set the minute .
	The setting flashes.
	Enter the minute with arrow buttons (0 to 59)
	and confirm with OK.

2x 	Back to Normal display.
--	-------------------------


6.2.5 Selecting the temperature unit

You can choose between degrees Celsius °C and degrees Fahrenheit °F.

If the temperature unit is changed, the temperature set-point and limits are converted accordingly.



Also when specifying the ramp function (see chap. 6.2.9) this setting is accordingly taken as the basis.

	C = degrees Celsius	0 °C = 31 °F	Conversion: [Value in °F] = [Value in °C] * 1.8 + 32
	F = degrees Fahrenheit	100 °C = 212 °F	

		From Normal display
without fan 5x with fan 6x		with the arrow-up button to the user menu
		Confirm with OK.
		Enter the password (factory setting: 00 00) and confirm each entry with OK.
		With the arrow-up button to the temperature unit selection menu.
		The current temperature unit is shown.
		Press OK to select the temperature unit.
		The setting flashes.
		Select the setting with arrow buttons
		and confirm with OK.
2x		Back to Normal display.

6.2.6 Set-point entry for temperature and fan speed

		From Normal display
		with the arrow-up button to the Set-point entry menu.
		The current temperature set-point is displayed.
		Press OK to enter the temperature set-point.
		The temperature set-point flashes.

	Enter the temperature set-point with arrow buttons with an accuracy of one degree (ED, FD, FED)
	and confirm with OK.
	Back to Normal display.
or	with chambers with adjustable fan speed (FED):
	go on to enter the fan speed.
	The fan speed set-point is displayed.
	Press OK to enter the fan speed
	The fan speed set-point flashes
	Adjust the fan speed in steps of 10 % with arrow buttons
	and confirm with OK.
	Back to Normal display.
	Check and/or adjust the safety controller following any changes of the set-point (chap. 6.2.9).

6.2.6.1 Set-point entry for temperature in two-door chambers (ED720, FED720)

With two-door chambers (ED720, FED720) the maximum adjustable temperature set-point depends on the ambient temperature:

- Range from 18 °C up to 26 °C: maximum temperature set-point: 300 °C
- Range from above 26 °C up to 40 °C: maximum temperature set-point decreasing with increasing ambient temperature


Figure: Maximum temperature set-point depending on the ambient temperature

Ambient temperature	Maximum temperature set-point
18 °C up to 26 °C	300 °C
28 °C	298 °C
30 °C	296 °C
32 °C	294 °C
34 °C	292 °C
36 °C	290 °C
38 °C	288 °C
40 °C	286 °C

This ensures the maximum lifetime of the controller.


6.2.7 Adjusting the air flap position

Opening the air flap in the exhaust duct serves to adjust the air change.


The position of the air flap in the exhaust duct serves to adjust the fresh air entry. With the open the air flap, fresh air can enter through the fresh air tube. For chambers with fan, the fan operation will increase fresh air entry.

If the air flap is completely open, the spatial temperature accuracy can be negatively influenced.

	From Normal display
without fan 2x with fan 3x	with the arrow-up button to the Adjusting the air flap position menu.
	The current air flap position is shown.
	Press OK to select the air flap position.

	The setting flashes.
	Select the position with arrow buttons
	and confirm with OK.
	Back to Normal display.

There are the following options:

Air flap closed	Air flap slightly opened	Air flap half opened	Air flap almost open	Air flap open
				
				

The set value can be set in 15° increments.

6.2.8 Changing the passwords for user level and general controller functions

In this menu you can change the passwords for access to the user menu and to all controller functions.

You can set two passwords for different access levels:

L1 (level 1): The password enables access control to the user level

L2 (level 2): The password enables access control to all controller functions

Factory setting for both passwords is 00 00 (no password assigned).

As soon as a password has been assigned, access to the respective functions is blocked and only available after entering the correct password.

	From Normal display
without fan 5x with fan 6x	 with the arrow-up button to the user menu
	Confirm with OK.
	Enter the password (factory setting: 00 00) and confirm each entry with OK.
2 x 	With the arrow-up button to the password setting menu .
	Password L1 for access to the user level. Confirm with OK.
	The current password level L1 is shown. The setting flashes. You can change between L1 and L2 with arrow buttons.
	Select the setting with arrow buttons (if desired)
	and confirm with OK.
	
	
	
	The current password for the selected password level is shown. The left two digits are flashing.
	Enter the desired numbers with arrow buttons,
	confirm with OK and go on.

	The right two digits of the password are flashing.
	Enter the desired numbers with arrow buttons
	and confirm with OK.
	The modified password (L1 or L2 depending on the selection) is shown (example: L1).

		If you want to change between L1 and L2 , confirm with OK. Thereafter you can change to the other password level and also modify the password.
2x		Back to Normal display.

	Keep in mind any modification of the password. There is no access to the user menu without the correct password L1. Without the correct password L2 access control to all controller functions is blocked.
---	--

6.2.9 Slow/ramp function

You can program temperature ramps in order to extend heating up times. This may be necessary in some cases to prevent temperature stress in the material during the heating up phase. Temperature ramps should only be used if required. Using them may result in considerably slowing down the heating up times. When the delayed (linear) heating function is disabled, the camera will carry out heating with maximum power.

The entry means the nominal value gradient and limits the maximum temperature increase to this value. Due to the heat and evaporation energy assumed by the drying material, smaller temperature gradients may also result.

A ramp proceeds from a previously entered set-point to a new, higher one. The temperature must be equilibrated to the start set-point.

Perform the setting in the following 3 steps:

1. Enter the temperature set-point as start slow/ramp set-point and let the temperature equilibrate to this value.
2. Define the temperature increase (slow/ramp gradient) in °C/min or in °F/min in the setting menu "Ramp function"

You can select a gradient from "0.0" up to "1.0" or from "1" up to "10" according to the oven type.

When setting the gradient to "0.0" or "0", the slow/ramp function is turned off. The chamber will then heat up with its maximum heating capacity.

000a will heat up in accordance with the entered gradient, i.e. with a speed of xx degrees per minute. A heating-up rate of 4 °C/min for the ED, FD and FED ovens can be regarded as a realistic maximum.


3. Enter the target slow/ramp set-point in the "ramp function" setting menu.

As soon as the entries have been adopted, the slow/ramp function is activated. The chamber heats up with the entered gradient, if the set ramp target value lies above the actual temperature value.

During ramp operation the effective slow/ramp set-point continually rises in accordance to the entered gradient from the previously entered set-point to the new one. The actual value follows this continually changing effective slow/ramp set-point. As soon as the heating-up target value is reached, this temperature is maintained constant.

The actual temperature value, the selected gradient, and the target value are shown in Normal display. The effective set value of delayed (linear) heating is a function of the set temperature value.


6.2.9.1 Setting and displaying the slow/ramp function

	From Normal display
without fan 3x with fan 4x 	with the arrow-up button to the Ramp function menu.
	Slow/ramp function (not programmed)
	Press OK to enter the gradient in degree per minute.
	The gradient flashes.
	Enter the gradient (from 0 to 9)
	Confirm with OK and go on to enter the target ramp set-point
	The target ramp set-point value flashes.
	Enter the target ramp set-point
	and confirm with OK.
	Selected ramp with gradient 1 and target ramp set-point 40 °C (example)
	Back to Normal display.
	Normal display showing the actual temperature, the gradient and the target ramp set-point (example: 40 °C)

Temperature course with slow/ramp function


 Normal display	 after 1 minute	 Normal display	 after 2 minutes	 Normal display	Etc.
---	--	---	--	---	------

6.2.9.2 Displaying the effective slow/ramp set-point and changing the target slow/ramp set-point

	From Normal display showing the actual temperature (example: 27 °C), gradient (example: 1) and the target ramp set-point (example: 40 °C)
	go to the display of the effective ramp set-point
	The effective ramp set-point is displayed
	Press OK to set the target ramp set-point
	The target ramp set-point value flashes.
	Enter the target ramp set-point
	and confirm with OK.
	The effective ramp set-point is displayed
2x 	Back to Normal display.
	Normal display showing the actual temperature, the gradient and the modified target ramp set-point (example: 50 °C)

6.2.9.2 Turning off the slow/ramp function

To turn off the ramp function, the gradient must be set to zero in the relevant menu. The set-point can be entered as desired.

	From Normal display
without fan 3x  with fan 4x	with the arrow-up button to the Ramp function menu
	Ramp function with programmed gradient and target ramp set-point
	Press OK to enter the gradient in degree per minute.





	The gradient flashes.
	Enter the gradient zero (for turning off the ramp function)
	Confirm with OK and go on to the ramp target value
	The target ramp set-point value flashes.
	Enter the target ramp set-point
	and confirm with OK.
	Slow/ramp function (not programmed)
	Back to Normal display.
	Normal display showing the actual temperature, the gradient and the target ramp set-point (example: 40 °C)
2x 	Back to Normal display.

Instead of turning off the ramp function (gradient = 0), you can also modify the settings of the gradient and target ramp set-point in this menu.

6.2.9.3 Turning off the slow/ramp function

To turn off the ramp function, the gradient must be set to zero in the relevant menu. The set-point can be entered as desired.

	From Normal display
without fan 3x  with fan 4x	with the arrow-up button to the Ramp function menu
	Ramp function with programmed gradient and target ramp set-point
	Press OK to enter the gradient in degree per minute.

	The gradient flashes.
	Enter the gradient zero (for turning off the ramp function)
	Confirm with OK and go on to the ramp target value
	The target ramp set-point value flashes.
	Enter the target ramp set-point

6.2.10 Timer functions

There are 3 timer functions:

	Timer function “Delayed Off” The selected timer countdown immediately starts running down. When the timer expires the heating turns off. Ovens with fan: The fan may continue working according to the selected settings. The maximum selectable time is 9:23:59 (Days: HH:MM)
	Timer function “Temperature dependent Delayed Off” The selected timer countdown only starts running down, when the actual value reaches or exceeds the selected set-point. When the timer expires the heating turns off. Ovens with fan: The fan may continue working according to the selected settings. The maximum selectable time is 9:23:59 (Days: HH:MM)
	Timer function “Delayed On” The selected timer countdown immediately starts running down, the heating turns off. When the timer expires the heating turns on and remains in continuous operation. Ovens with fan: The fan may be working according to the selected settings. The maximum selectable time is 9:23:59 (Days: HH:MM)

FED ovens offer all three timer functions.

ED and FD ovens offer the timer function “Delayed Off”. The maximum selectable time is 23:59 (HH:MM)

Stage	Heating	Fan (ovens with fixed fan speed: FD)	Fan (ovens with adjustable fan speed: FED)
Timer function "Delayed Off"			
Timer running	Control to the temperature set-point	On (100 %)	Rotation speed according to fan speed set-point
After the timer expired	Off	On (100 %) or Off (0 %) acc. to selection	Rotation speed according to setting of timer functions
Timer function "Temperature dependent Delayed Off"			
Possibly heating-up phase until the temperature set-point is reached	Control to the temperature set-point	---	Rotation speed according to fan speed set-point
Timer running	Control to the temperature set-point	---	Rotation speed according to fan speed set-point
After the timer expired	Off	---	Rotation speed according to setting of timer functions
Timer function "Delayed On"			
Timer running	Off	---	Rotation speed according to setting of timer functions
After the timer expired	Control to the temperature set-point	---	Rotation speed according to fan speed set-point

General information on the setting:


In the setting menus of the timer functions, it is always required to confirm all parameters with OK, otherwise all entries made will be lost.

Timer run-time is set in days, hours, and minutes. If days have been entered, they are shown in the controller display preceded by an underscore:

	Setting: 0 days (not shown), 10 hours, 30 minutes		Setting: 2 days, 10 hours (minutes not shown)
---	---	--	---

6.2.10.1 Timer function "Delayed Off"


6.2.10.1.1 Entry and activation of the timer run-time and fan setting

	From Normal display
	with the arrow-down button get to the Timer functions “Delayed Off” menu (with connected USB device: press the arrow-down button twice)
	Current Timer function “Delayed Off”
	Confirm with OK and go on to enter the days of the timer run-time .
	The current timer run-time (days) is shown. The day value flashes.
	Enter the days of the timer run-time (from 0 to 9)
	Confirm with OK and go on to enter the hours of the timer run-time .
	The hour value flashes.
	Enter the hours of the timer run-time (from 0 to 23)
	Confirm with OK and go on to enter the minutes of the timer run-time
	The minute value flashes.
	Enter the minutes of the timer run-time (from 0 to 59)
	and confirm with OK.

Ovens without fan (ED):

	Timer function “Delayed Off” is activated
	Back to Normal display.

Ovens with fixed fan speed (FD): Select whether the fan shall operate or not after the timer expired	
	The current setting of the fan operation is shown
	Select fan operation: On (100 %) or Off (0 %)
	and confirm with OK.
	Timer function "Delayed Off" is activated
	Back to Normal display.

Ovens with adjustable fan speed (FED): Enter the fan speed set-point valid for the time after the timer expired.	
	The current fan speed set-point is shown.
	Enter the fan speed in steps of 10 %: 0 % and 40 % up to 100 %
	and confirm with OK.
	Timer function "Delayed Off" is activated
	Back to Normal display.

Normal display during timer operation with timer function "Delayed Off"

	The timer run-time is running until turning off the heating.
---	--

Normal display after the timer expired:

	Ovens without fan (ED): Timer function "Delayed Off". The timer has expired. The heating is off.
	Ovens with fixed fan speed (FD): Timer function "Delayed Off". The timer has expired. The heating is off. The fan is operating ("On") or not ("Off", вкл), as selected
	Ovens with adjustable fan speed (FED): Timer function "Delayed Off" (Отключение с задержкой по времени). The timer has expired. The heating is off. The fan is operating at the selected fan speed.

When the timer has expired the heating is off. The oven cools down to ambient temperature.

To restart the oven, you need to turn off the timer function (Chapter 6.2.10.1.2)

6.2.10.1.2 Turning off the timer functions or changing the settings

To turn off the timer function "Delayed Off" during the time when it is still running, set all time values (days, hours, minutes) to zero in the relevant menu (terminating the timer function). As long as the timer is running, the timer function settings can be subsequently modified in this menu.

After the timer expired you can turn off the timer function by pressing the OK button. Alternatively, you can also set the time to zero in the relevant menu.

Turning off the timer function when the timer has expired

In Normal display press the OK button.

Alternatively, you can set the time to zero:

	From Normal display
	with the arrow-down button get to the Timer functions "Delayed Off" menu (with connected USB device: press the arrow-down button twice)
	Ovens without fan (ED): Timer function "Delayed Off", The timer has expired. The heating is off.
или	
	Ovens with fixed fan speed (FD): Timer function "Delayed Off", The timer has expired. The heating is off. The fan is operating ("On") or not ("Off"), as selected
или	
	Ovens with adjustable fan speed (FED): Timer function "Delayed Off", The timer has expired. The heating is off. The fan is operating at the selected fan speed.

	Confirm the timer function “Delayed Off” with OK.
see Chapter 6.2.10.1.1	Set all time values to zero
2x 	Back to Normal display.








Changing or terminating the timer function when the timer is running:

	From Normal display
	With the arrow-down button get to the Timer functions “Delayed Off” menu (with connected USB device: press the arrow-down button twice)
	Ovens without fan (ED): Timer function “Delayed Off”, the timer is running.
или	
	Ovens with fixed fan speed (FD): Timer function “Delayed Off”, the timer is running.
или	
	Ovens with adjustable fan speed (FED): Timer function “Delayed Off”, the timer is running.
	Confirm the timer function “Delayed Off” with OK.
see Chapter 6.2.10.1.1	Modify the time values or set them to zero to terminate the timer functions
2x 	Back to Normal display.

6.2.10.2 Timer function “Temperature dependent Delayed Off” (FED ovens)

6.2.10.2.1 Entry and activation of the timer run-time, fan setting and set-point entry

	From Normal display
3 x 	with the arrow-down button get to the Timer functions “Temperature dependent Delayed Off” menu (with connected USB device: press the arrow-down button 4 times)
	Current Timer function “Temperature dependent Delayed Off”
	Confirm with OK and go on to enter the days of the timer run-time .

	The current timer run-time (days) is shown. The day value flashes.
	Enter the days of the timer run-time (from 0 to 9)
	Confirm with OK and go on to enter the hours of the timer run-time .
	The hour value flashes.
	Enter the hours of the timer run-time (from 0 to 23)
	Confirm with OK and go on to enter the minutes of the timer run-time
	The minute value flashes.
	Enter the minutes of the timer run-time (from 0 to 59)
	Confirm with OK and go on to enter the temperature set-point
	The temperature set-point flashes. When the set-point is reached, the timer starts running down.
	Enter the temperature set-point,
	Confirm with OK and go on to enter the fan speed set-point valid for the time after the timer has expired
	The current fan speed is shown
	Enter the fan speed in steps of 10 %: 0 % and 40 % up to 100 %
	and confirm with OK.
	Timer function "Temperature dependent Delayed Off" is activated
	Back to Normal display.

Normal display during timer operation with timer function “Temperature dependent Delayed Off”

	Normal display while the timer is running If the entered set-point exceeds the current actual value, the chamber heats up.
	
	The defined time only begins to run when the current value is by 1 °C below the set point. After the defined time has expired, the heating is turned off. Normal display during the timer run-time. The timer run-time until turning off the heating (and possibly the fan) is running.

Normal display after the timer expired:

	Timer function “Temperature dependent Delayed Off” The timer has expired. The heating is off. The fan is operating at the selected fan speed.
---	---

When the timer has expired the heating is off. The chamber cools down to ambient temperature.

To restart the oven, you need to turn off the timer function (Chapter 6.2.10.2.2).

6.2.10.2.2 Turning off the timer functions or changing the settings

To turn off the timer functions “Temperature dependent Delayed Off” during the time when it is still running, set all time values (days, hours, minutes) to zero in the relevant menu (terminating the timer function). As long as the timer is running, the timer function settings can be subsequently modified in this menu.

After the timer expired you can turn off the timer function by pressing the OK button. Alternatively, you can also set the time to zero in the relevant menu.

Turning off the timer functions when the timer has expired

	From Normal display
3 x 	with the arrow-down button get to the Timer functions “Temperature dependent Delayed Off” menu (with connected USB device; press the arrow-down button 4 times)
	Timer function “Temperature dependent Delayed Off”, The timer has expired
	Confirm the timer function “Temperature dependent Delayed Off” with OK.
see Chapter 6.2.10.2.1	Set all time values to zero
	Back to Normal display.

Changing or terminating the timer function while the timer is running:



		From Normal display
3 x		with the arrow-down button get to the Timer functions “Temperature dependent Delayed Off” menu (with connected USB device: press the arrow-down button 4 times)
		Timer function “Temperature dependent Delayed Off”, the timer is running
		Confirm the timer function “Temperature dependent Delayed Off” with OK.
see Chapter 6.2.10.2.1		Modify the time values or set all time values to zero to terminate the timer functions
		Back to Normal display.

6.2.10.3 Timer function “Delayed On” (FED ovens)

6.2.10.3.1 Entry and activation of the timer run-time and fan setting

		From Normal display
2 x		with the arrow-down button get to the Timer functions “Delayed On” menu (with connected USB device: press the arrow-down button 3 times)
		Current Timer function “Delayed On”
		Confirm with OK and go on to enter the fan speed set-point valid during the time the timer is running
		The current fan speed set-point is shown.
		Enter the fan speed in steps of 10 %: 0 % and 40 % up to 100 %
		Confirm with OK and go on to enter the days of the timer run-time.
		The current timer run-time (days) is shown. The day value flashes.
		Enter the days of the timer run-time (from 0 to 9)
		Confirm with OK and go on to enter the hours of the timer run-time.

	The hour value flashes.
	Enter the hours of the timer run-time (from 0 to 23)
	Confirm with OK and go on to enter the minutes of the timer run-time

	The minute value flashes.
	Enter the minutes of the timer run-time (from 0 to 59)
	Confirm with OK and go on to enter the temperature set-point
	The temperature set-point flashes.
	Enter the temperature set-point,
	Confirm with OK and go on to enter the fan speed set-point valid for the time after the timer has expired
	The current fan speed set-point is shown.
	Enter the fan speed in steps of 10 % - 0 % and 40 % up to 100 %
	and confirm with OK.
	Timer function "Delayed On" is activated
2x 	Back to Normal display.

Normal display during timer operation with timer function "Delayed On"

	The Timer run-time until turning on the heating is running. Timer function "Delayed On". The heating is off, temperature approximates ambient temperature.
---	--

Normal display After the timer expired:

	The timer has expired. The timer function is off. The heating is active to equilibrate the temperature set-point. The fan is operating at the selected fan speed.
---	--

6.2.10.3.2 Changing the settings

After the timer expired, the timer function “Delayed On” deactivates, therefore turning the function off is not required.

As long as the timer is still running, the timer function settings can be subsequently modified in this menu.

To terminate the timer function, all time values (days, hours, minutes) must be set to zero in the relevant menu.

Changing or terminating the timer function when the timer is running:

	From Normal display during timer operation
2 x 	with the arrow-down button get to the Timer functions “Delayed On” menu (with connected USB device: press the arrow-down button 3 times)
	Timer function “Delayed On”, the timer is running. The fan is operating at the selected fan speed.
	Confirm the timer function “Delayed On” with OK.
see Chapter 6.2.10.3.1	Modify the time values or set all time values to zero to terminate the timer function
2x 	Back to Normal display.

6.2.10.4 Temperature programming example (FED ovens)

The oven shall heat up to a temperature of 50 °C, maintain this temperature for three hours and then turn off.

Proceeding: Select the timer function “Temperature dependent Delayed Off “(Chapter 6.2.10.2) and perform the following settings:

- Enter a timer run-time of 3 hours
- Enter the set point: 50 °C
- Specify the fan speed after the timer expires

7. Overtemperature protection (independent thermal protection system)

An independent thermal protection system is a built-in element that interrupts the operation of the oven in case of overheating and protects the oven, its surroundings and the load from unacceptable excessive temperatures.

It is recommended to check the functionality of the independent thermal protection system at appropriate intervals. It is recommended to carry out such a check, for example, before starting a long working cycle, and it is necessary to involve personnel specially authorized to carry out this type of work.

Please follow the safety rules and internal regulations in the laboratory and medical facilities.

Depending on the model, in accordance with DIN 12880:2007, cabinets are equipped with an independent thermal protection system of class 1 (overheating device), an independent thermal protection system of class 2 (temperature limiter), and also, upon customer request, an independent thermal protection system of class 3.1 (temperature protection).

7.1 Independent thermal protection system for E28 ovens

E28 ovens are equipped with an overheating protection device acc. with DIN 12880:2007 (independent thermal protection system class 1), which serves the purpose of protecting the cabinet from damage that may be caused by overheating.

If the cut-off temperature 320 °C is reached the thermal protection device will turn off the equipment

The safety shut-off device is located inside the device. Only a service technician can replace it. Therefore, please contact an authorized service center or the BINDER service department.

7.2 Independent thermal protection system for ED23, FD23, ED400, FED400 ovens

7.2.1 Independent thermal protection system Class 2

The Independent Thermal Protection System Class 2 is a safety device (electronic temperature limiter with mechanical reset). Upon request by the consumer and only in the factory, the independent thermal protection system of class 2 can be replaced by an independent thermal protection system of class 3.1.

Operating principle:

The independent thermal protection system class 2 is functionally and electrically independent of the cabinet management system (controller) and shuts down the cabinet completely at all terminals. If a problem occurs in the temperature protection device, the independent thermal protection system disables the cabinet control system (controller). This status is indicated by the indicator (7a), and if installed optionally, an acoustic alarm is also activated using an audible signal. The operation of the independent thermal protection system is checked by slowly turning the control knob (7) counterclockwise until the

system is turned off. The independent thermal protection system is switched back on by pressing the reset button (7b).


Figure: Class 2 independent thermal protection system

Functioning:

When the control knob (7) is set to the end position, the class 2 independent thermal protection system acts as an equipment protection device. If it is set to a temperature higher than the temperature selected on the controller, it acts as a safety device for the material being processed.



When the independent thermal protection system is triggered, the red signal lamp (7a) starts flashing, and with an optional setting, an acoustic signaling is additionally triggered using an acoustic signal. To restart the oven, follow these steps:

- Disconnect the oven from the power supply
- Determine and eliminate the cause of the malfunction
- Start the system of independent thermal protection class 2 by pressing the reset button (7b).
- Restart the oven as described in chapter 5.2.1.

Setting:

To check to what temperature the system of independent thermal protection class 2 is triggered, turn on the device and set the desired fixed value on the controller.

The scale divisions 1 to 10 correspond to a temperature range of 30°C to 320°C and are intended to aid installation.

1. Set the control knob (7) of the independent thermal protection system to its end position (position 10). Use a coin to install.	
2. When the fixed temperature is reached, move the control knob (7) to the trigger point (move counterclockwise).	
3. The trigger point is identified by the red warning lamp (7a) lighting up; the reset button (7b) is pushed forward.	
In the presence of an acoustic signaling and when the acoustic signal is on (Chapter 12.6), the buzzer sounds as an additional signal. It can be turned off by switch (10) or (9) (Fig. 6, item 2.2.2).	
4. The optimal setting of the independent thermal protection system is achieved by moving the control knob clockwise in a circle by one scale division.	
5. Press the reset button (7b) again.	

The independent thermal protection system is only active when the reset button (7b) is pressed.

When the class 2 independent thermal protection system is triggered, the red signal lamp (7a) comes on, the reset button (7b) is pushed out and the device switches off reliably in all terminals.

The settings of the class 2 independent thermal protection system should be regularly checked and, if necessary, adjusted with each change in the setpoint temperature.

7.2.2 Independent thermal protection System class 3.1

The independent thermal protection system of class 3.1 is a protective device (electronic temperature protection with mechanical restart) and is a cabinet component that is installed at the customer's request in the factory instead of an independent thermal protection system of class 2.

Operating principle:

The Class 3.1 independent thermal protection system is a functionally and electrically independent cabinet control system (controller); in the event of an error, an adjustment function can be used.


Figure: Class 3.1 independent thermal protection system

When the control knob is set to the end position, the class 3.1 independent thermal protection system acts as an equipment protection device. If it is set to a temperature higher than the temperature selected on the controller, it acts as a safety device for the material being processed.


When the independent thermal protection system is triggered, the red signal lamp (7a) starts flashing, and with an optional setting, an acoustic signaling (option, Chapter 12.6) is additionally triggered using an acoustic signal. To restart the oven, follow these steps:

- Disconnect the oven from the power supply
- Determine and eliminate the cause of the malfunction
- Restart the oven (chapter 5.2.1.).

Setting:

To check to what temperature the system of independent thermal protection class 3.1 is triggered, turn on the oven and set the required temperature setpoint on the controller.

The scale divisions 1 to 10 correspond to a temperature range of 63°C to 350°C and are intended to aid installation.

1. Turn the control knob (7) of the system of independent thermal protection class 3.1 with a coin to the end position (equipment protection).	
2. When the setpoint is reached, turn the control knob (7) to the shutdown point (turn it counterclockwise).	
3. The switch-off point is identified by the lighting up of the red warning light (7a). If the optional acoustic signaling is present and the buzzer is activated (chap. 12.6), the buzzer gives an additional signal. It can be turned off by switch (10) or (9) (Fig. 6, item 2.2.2).	
4. To achieve the optimal setting of the system of independent thermal protection class 3.1, turn the control knob clockwise by about one scale division, while the red signal lamp (7a) goes out.	

	The settings of the class 3.1 independent thermal protection system should be regularly checked and, if necessary, adjusted with each change in the setpoint temperature.
--	--

7.3 Independent thermal protection system for ovens ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720

7.3.1 Independent thermal protection system Class 1

As standard, the cabinets are equipped with an overheating protection device acc. with DIN 12880:2007 (independent thermal protection system class 1). It serves the purpose of protecting the oven, its installation location and its contents from exceeding the maximum permissible temperature values. If the temperature exceeds 20°C - 30°C above the nominal temperature of the chamber, the message "Overtemperature" will appear on the controller display and the system of independent thermal protection class 1 turns off the heating.

Shutdown temperature:

For ED 240, ED 720: 320°C

For ED 53, ED 115, FD 53, FED 53: 330°C

For FD 115, FD 240, FED 115, FED 240, FED 720: 350°C

To restart the oven, follow these steps:

Disconnect the oven from the power supply.

Determine and eliminate the cause of the malfunction.

Let the oven cool down.

Restart the oven.

As soon as the temperature inside the chamber after restarting falls below the response temperature of the overtemperature protection class 1, the alarm will be cleared automatically.

Restart temperature:

For ED 240, ED 720: 220°C

For ED 53, ED 115, FD 53, FED 53: 230°C

For FD 115, FD 240, FED 115, FED 240, FED 720: 250 °C

7.3.2 Independent thermal protection system class 2 ("temperature limiter")

The Independent Thermal Protection System Class 2 is a safety device (electronic temperature limiter with electronic reset). Upon request by the consumer and only in the factory, the independent thermal protection system of class 2 can be replaced by an independent thermal protection system of class 3.1.

An independent thermal protection system of class 2 is used to limit the temperature selected by the consumer inside the cabinet from the set value on the independent thermal protection system device.

The limit temperature is entered as the set value of the independent thermal protection system, i.e. as the absolute maximum allowable temperature.

If the maximum temperature is exceeded, the independent thermal protection system completely turns off the heating until it is manually restarted. This is indicated visually by warning messages and acoustic signaling using an audible signal when the buzzer is activated.

If the system of independent thermal protection class 2 has turned off the heating, it is recommended to proceed as follows:

- Disconnect the oven from the power supply.
- Determine and eliminate the cause of the malfunction.
- Restart the oven.
- Clear alarm message.

7.3.3 Independent thermal protection system class 3.1 ("temperature protection" - protection against overtemperature)

The independent thermal protection system of class 3.1 is a protective device (electronic temperature protection with automatic restart) that is included in the cabinet, which is installed at the customer's request in the factory instead of an independent thermal protection system of class 2.

An independent thermal protection system of class 3.1 is used to protect against the temperature selected by the consumer exceeding the temperature inside the cabinet from the set value on the independent thermal protection system device.

The limit temperature is entered as the set value of the independent thermal protection system, i.e. as the absolute maximum allowable temperature.

If the maximum temperature is exceeded, the independent thermal protection system completely turns off the heating until it is manually restarted. This is indicated visually by warning messages and acoustic signaling using an audible signal when the buzzer is activated.

The independent thermal protection system ensures that the oven is controlled until the selected chamber temperature cools below the set value on the device of the independent thermal protection system.

If the system of independent thermal protection class 3.1 has taken over the control of the oven, it is recommended to proceed as follows:

- Disconnect the oven from the power supply.
- Determine and eliminate the cause of the malfunction.
- Restart the oven.
- Clear alarm message

7.3.4 Setting the safety controller set-point

A limit temperature is entered as the safety controller set-point, i.e. the absolute maximum permitted temperature value.

Example: Temperature set-point 45 °C, safety controller set-point 50 °C.

	Regularly check the safety controller setting relating to the entered temperature set-point. Set the safety controller set-point by approx. 2 °C to 5 °C above the desired temperature set-point.
---	--

	From Normal display
without fan 4x with fan 5x	 with the arrow-up button to the Safety controller set-point setting menu.
	The current safety controller set-point is shown (class 2 "temperature limiter" or class 3.1 "temperature protection" depending on the chamber type).
	Press OK to enter the safety controller set-point.
	The safety controller set-point flashes.
	Enter the safety controller set-point with arrow buttons: 10 °C up to 300 °C (with an accuracy of one degree) with ED, FD, FED
	and confirm with OK.
	Back to Normal display.

7.3.5 Alarm message and proceeding in case of an alarm

The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 2.3.6) the buzzer sounds.

Safety controller class 2 ("temperature limiter")

The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 2.3.6) the buzzer sounds.

The heating turns off.

Resetting the alarm:

With option audible alarm with activated buzzer: Mute the buzzer pressing the OK button.

As soon as the inner chamber temperature has cooled down below the safety controller set-point, the alarm icon is lit permanently. You can reset the alarm message on the controller. To do this, reset the alarm message in the safety controller menu with the OK button. The heating is released and temperature control is resumed by the controller.

Safety controller class 3.1 ("temperature protection")

The alarm icon flashes in the display. In case of the option audible alarm with activated buzzer (chap. 2.3.6.) the buzzer sounds.

The heating turns off.

Resetting the alarm:

With option audible alarm with activated buzzer: Mute the buzzer pressing the OK button.

As soon as the inner chamber temperature has cooled down below the safety controller set-point, you can reset the alarm message in the safety controller menu with the OK button. The heating is released and temperature control is resumed by the controller.

Note:

When the safety controller class 2 or class 3.1 had been activated, you should disconnect the chamber from the power supply and have an expert examine and rectify the cause of the fault.

8. Data recording via USB interface for ED53, ED115, ED240, ED720, FD53, FD115, FD240, FED53, FED115, FED240, FED720 ovens

The USB interface located in the triangular instrument panel serves to read out the measured values, which are put out in real time. Following data is recorded: Timer, actual temperature value, temperature set-point, object temperature senso (for ovens having Object temperature display as an option), fan speed (for ovens with fan), air flap position, safety controller, analog output (option), heating regulation ratio.

Connect only USB sticks to the USB interface.

Data are stored directly on the USB stick. They are issued in the selected language as a spreadsheet with the file extension ".csv" and can be further processed in the desired program.

8.1 Starting data recording


Connect the USB stick to the interface located in the triangular instrument panel.

	From Normal display
	with the arrow-down button get to the USB menu
	Current state of data recording: stopped
	Confirm with OK.

Chambers without the real time clock option do not save date and time after shut-down. To make sure that recorded data is provided with the correct date and time, with these chambers it is required to enter date and time first:

	The current date is displayed. For setting, proceed as described in chap. 6.2.4.
---	---

After entry of the minute and confirmation with OK, further setting in the USB menu continues.

	The data recording interval flashes.
	Enter the interval in minutes
	and confirm with OK.
	The current state of data recording (stopped) flashes.
	Change to start data recording.
	The current state of data recording (started) flashes.
	Confirm with OK.
	Data recording is running.
	Back to Normal display.

Data recording continues also during stand-by mode of the chamber. Disconnecting the power supply interrupts data recording, which continues after the power returns. To terminate data recording, stop it via the menu (chap.8.2).

8.2 Terminating data recording

	From Normal display
	with the arrow-down button get to the USB menu
	Current state of data recording: running
	Confirm with OK.
	The current state of data recording (running) flashes
	Change to stop data recording
	The current state of data recording (stopped) flashes
	Confirm with OK.
	Data recording is stopped.
	Back to Normal display.

9. Network configuration for FED53, FED115, FED240, FED720 ovens with Ethernet interface

The settings of this submenu are used to network cabinets via an Ethernet interface, for example, to connect them to the APT-COM communication software supplied by BINDER.

This menu subsequently offers the following settings:

- Display of the chamber's **MAC address** (no setting)
- Switching on and off the **DHCP state**
- Entering the **IP address**
- Entering the **subnet** mask number
- Entering the **gateway** number

		From Normal display
without fan 5x with fan 6x		with the arrow-up button to the user menu
		Confirm with OK.
		Enter the password (factory setting: 00 00) and confirm each entry with OK.
5 x		with the arrow-up button to the Ethernet menu.
		Confirm with OK.
		The first digits of the MAC address are shown.
		Continue with OK. The middle digits of the MAC address are shown.
		Continue with OK. The last digits of the MAC address are shown.

Displaying the chamber's MAC address serves to identify the chamber in the Ethernet network.

Example: **00-0F-67-0F-42-40**

00-0F	67-0F	42-40
		

	Continue with OK to switching on and off the DHCP state
	The current DHCP state is displayed (on). The setting flashes.
	Select the setting with the arrow buttons.

There are the following options:

DHCP switched on:	DHCP switched off:
	

To configure the network settings, the DHCP state must be switched off. Otherwise, the DHCP-server would assign the network configuration..

	Confirm with OK and go on to enter the IP address .
	The first digits of the current IP address are shown. The setting flashes.
	Enter the desired values with the arrow buttons
	and confirm with OK. The next digits of the IP address are shown. Enter the desired values accordingly.

Example value: 0.0.0.0

0	0	0	0
			

	Confirm with OK and go on to enter the subnet mask .
	The first digits of the current subnet mask are shown. The setting flashes.
	Enter the desired values with the arrow buttons
	and confirm with OK. The next digits of the subnet mask are shown. Enter the desired values accordingly.

Example value: 255.255.255.0

255	255	255	0
			

	Confirm with OK and go on to enter the gateway .
	The first digits of the current gateway are shown. The setting flashes.
	Enter the desired values with the arrow buttons
	and confirm with OK. The next digits of the gateway are shown. Enter the desired values accordingly.

Example value: 0.0.0.0

0	0	0	0
			

	Confirm with OK.
2x 	Back to Normal display.

10. Cleaning and decontamination

Clean the chamber after each use to avoid possible corrosion damage by ingredients of the loading material.

Prior to renewed startup, allow the chamber to completely dry after all cleaning and decontamination measures.

	 DANGER
	Danger of electric shock. Danger to life. - Ø DO NOT pour water or detergent onto interior or exterior surfaces. - Ø Do NOT put ANY cleaning aids (cloth or brush) into slots or openings on the chamber ➤ Unplug this product before cleaning. Pull out the power connector. ➤ Dry the device completely before turning it on.

10.1 Cleaning oven

In case of contamination of the chamber with biologically or chemically hazardous substances, it is necessary to perform proper disinfection of the oven. Before disinfection, it is necessary to disconnect the oven from the mains by removing the power connector from the outlet. Do not use disinfectants that may cause harm due to reaction with the components of the device and the test material. Follow the operating

instructions and safety recommendations. Periodically check the internal condition of the camera. Disinfection is carried out as contamination occurs or at least once a week. Use a spray for chemical disinfection.

Recommended disinfectants:

Exterior surfaces Inner chamber Racks Door gasket	Use cleaning products with a neutral pH of 6.0-8.0, which do not destroy the structure of materials, including solutions containing at least 60-80% isopropyl or ethyl alcohol. Do not use cleaning products containing acids and halides. It is recommended to use a neutral detergent - "deconex 36 INTENSIVE" (composition: fatty alcohol ethoxylate ≥ 15 – $<30\%$, nonionic surfactants ≥ 15 - $<30\%$, phosphones $<5\%$), or another product with a similar composition.
Controller panel	Use cleaning products with a neutral pH of 6.0-8.0, which do not destroy the structure of materials. Do not use cleaning products containing acids and halides. It is recommended to use a neutral detergent - "deconex 36 INTENSIVE" (composition: fatty alcohol ethoxylate ≥ 15 – $<30\%$, nonionic surfactants ≥ 15 - $<30\%$, phosphones $<5\%$), or another product with a similar composition.
Zinc coated hinge parts, rear housing wall	It is not recommended to use a neutral detergent "deconex 36 INTENSIVE" or another product with a similar composition on galvanized surfaces.

	Do not use cleaning agents that may be harmful due to reactions with the components of the device and the material under test. If in doubt about the suitability of the detergent, please contact BINDER Service.
---	---

	Any corrosive damage that may occur due to the use of other cleaning agents is excluded from the BINDER GmbH warranty. Any corrosive damage caused by lack of cleaning is excluded from the responsibility of BINDER GmbH.
---	---

	Chlorides may be present in the foam and therefore cannot be used for cleaning.
---	---

	NOTICE Danger of corrosion. Risk of damage to the device. - Ø Do not use acid or chlorine reagents. - Ø Do not use neutral cleaning agents on some surfaces, eg galvanized hinge parts, back of the housing).
---	---

	In case of cleaning take care of adequate personal safety measures.
---	---

	Recommended safety measures: ➤ Wear tight-fitting safety goggles to protect your eyes. Suitable protective gloves for full contact: butyl rubber or nitrile rubber, burst time: > 480 minutes ➤ Neutral cleaning agent may be harmful to health if it comes into contact with the skin or is swallowed. Observe the instructions for use and safety recommendations written on the neutral detergent vials.
---	---

	 CAUTION Skin contact, ingestion. Damage to skin and eyes due to chemical burns. - Ø DO NOT swallow. Avoid contact with food and drink. - Ø Do not empty into drains. - Ø Wear protective gloves and goggles. - Ø Avoid skin contact with the neutral cleaning agent.
---	---

10.2 Decontamination

In case of contamination of the chamber with biologically or chemically hazardous substances, it is necessary to perform proper disinfection of the oven. Before disinfection, it is necessary to disconnect the oven from the mains by removing the power connector from the outlet. Do not use disinfectants that may cause harm due to reaction with the components of the device and the test material. Follow the operating instructions and safety recommendations. Periodically check the internal condition of the camera. Disinfection is carried out as contamination occurs or at least once a week. Use a spray for chemical disinfection.

Recommended disinfectants:

Inner chamber	Use cleaning products with a neutral pH of 6.0-8.0, which does not destroy the structure of materials, including solutions containing at least 60-80% isopropyl or ethyl alcohol or 0.25% copper sulfate solution (copper sulfate). Do not use cleaning products containing acids and halides It is recommended to use deconex SOLARSEPT disinfectant spray (composition: Propane-1-ol (also propanol, n-propyl alcohol) – <30%, Propane-2-ol (Isopropyl alcohol) - ≥15–<30%, N-(3-Aminopropyl)-N-dodecylpropane-1,3-diamine - <1% is well soluble in alcohol and water), or another product with a similar composition.
---------------	--

In case of interior contamination with hazardous biological or chemical substances, there are two possible procedures depending on the type of contamination and loading material:

1. Sterilization at 180°C for at least 30 minutes. Any flammable materials must be removed from the inner chamber before the sterilization.
2. Spraying (spray, spray gun) in the inner chamber with recommended disinfectants.

Before start-up, the device must be absolutely dry and ventilated, because explosive gases may form during the decontamination process.

3. If case of heavy contamination, the inner chamber parts must be removed (by a qualified service technician) for cleaning or replacement. After that, sterilize the inner chamber parts in a sterilizer or an autoclave. Racks can also be removed and sterilized.


	After using the spray, dry and ventilate the chamber.
---	---

	Recommended safety measures: ➤ Wear tight-fitting safety glasses to protect your eyes. Disinfectant spray solution may cause eye damage due to chemical burns. ➤ Follow the instructions for use and safety recommendations written on the disinfectant spray bottles.
---	--

 	 CAUTION Eye contact. Eye damage due to chemical burns. - Ø Do not empty into drains. - Ø Wear safety glasses.
---	--

11. Maintenance, troubleshooting

11.1 Maintenance intervals

	 DANGER Electrical hazard during live maintenance work. Deadly electric shock. - Ø The chamber must NOT become wet during operation or maintenance works. - Ø Do NOT remove the rear panel of the chamber. ➤ Disconnect the chamber before conducting maintenance work. Turn off the main power switch and pull the power plug. ➤ Make sure that all maintenance work will be conducted by licensed electricians or experts authorized by BINDER.
---	---

Ensure regular maintenance work is performed at least once a year.

	The warranty becomes void if maintenance work is conducted by non-authorized personnel.
---	---

	Replacement of the door seal is possible only in a cooled chamber. Otherwise, the door seal will be damaged.
---	--

We recommend taking out a maintenance agreement. Please consult BINDER Service.

BINDER telephone hotline:	+49 (0) 7462 2005 555
BINDER fax:	+49 (0) 7462 2005 93555
BINDER e-mail:	service@binder-world.com
BINDER service hotline Russia and CIS	+7 495 988 15 16
BINDER Internet website	http://www.binder-world.com
BINDER address	BINDER GmbH, Im Mittleren Oesch 5, 78532 Tuttlingen

International customers, please contact your local BINDER distributor.

11.2 Troubleshooting

Fault description	Possible cause	Required measures
Temperature		
Set-point temperature is not reached after specified time.	Door gasket defective.	Replace door gasket,
	Controller not adjusted.	Calibrate and adjust controller.
	Wrong voltage.	Check the power supply for correct voltage (chap. 4.2.)
<i>For FD, FED:</i> The fan doesn't turn or turns too slowly.	<i>For FED:</i> Fan speed set too low	<i>For FED:</i> Set fan speed to 100%
	Fan defective.	Contact BINDER service.
Chamber heating permanently, set-point not held.	Controller defective.	Contact BINDER service.
	Semiconductor relay defective	
	Controller not adjusted.	Calibrate and adjust controller.
Chamber doesn't heat up. Heating icon is displayed or red heating indicator light is on	Heating element defective.	Contact BINDER service.
	Semiconductor relay defective.	
Chamber doesn't heat up. No heating icon in the display. Controller display working. <i>For ED23, FD23, ED400, FED400:</i> red heating indicator is not shown	Timer run off.	Re-program the timer or turn it off. <i>For ED23, FD23, ED400:</i> Change its mode to "Continuous operation" <i>For FED400:</i> Set a new timer set-point
	Semiconductor relay defective.	Contact BINDER service.
	Controller defective.	
Chamber without function, the standby icon is displayed <i>For ED23, FD23, ED400, FED400:</i> red indicator of "Standby" mode is shown	Chamber in standby mode.	Hold the standby button down until the display lights up. <i>For ED23, FD23, ED400, FED400:</i> Hold the power on/off button (5) down until the display lights up.
<i>For ED, FD, FED:</i> Chamber without function. Alarm message "TLim" is displayed <i>For ED23, FD23, ED400, FED400:</i> The red signal lamp of independent thermal protection is on (7a).	Safety device class 2 has turned off the chamber.	Let the chamber cool down. Check the settings of the temperature set-point and independent thermal protection of the safety device class 2, press the RESET button (if any). If appropriate, select suitable limit value.
	Safety device class 2 defective.	Contact BINDER service.
<i>For ED23, FD23, ED400, FED400:</i> The inside temperature is too high. The red signal lamp of independent thermal protection is on (7a).	Independent thermal protection system, class 3.1 (optional) tripped.	Check the set points of temperature and of class 3.1. unit

Chamber without any function.	No power supply.	Check connection to power supply.
	For ED23, FD23, ED400, FED400: The fuse is blown.	For ED23, FD23, ED400, FED400: Check the fuse and replace it if necessary. If it blows again, contact BINDER Service.
	Controller defective.	Contact BINDER service.
Deviations from the indicated heating-up times.	Chamber fully loaded.	Load the chamber less or consider longer heating-up times.
Controller		
Message „1999” in the controller display	Sensor rupture between sensor and controller.	Contact BINDER service.
For ED23, FD23, ED400, FED400: The controller returns to Normal Display from any level	Buttons have not been pressed for more than approx. 30 seconds.	Repeat entering values, enter values faster.
Total		
For E28: The camera is not working The green indicator light (3) is off. The timer is set ∞ or for the selected time	There is no power supply	Check the connection to the power supply Check if the camera is turned on using the main power switch or timer.
	Incorrect voltage	Check that the voltage in the power supply is correct.
	The camera's fuse went off.	Check the camera fuse and replace it if necessary. If it works again, contact the BINDER customer service.
	The timer has expired	Check the timer setting
	The camera is defective	Contact the BINDER customer service.

Repairs should only be carried out by qualified personnel authorized by BINDER. Repaired equipment must comply with the requirements described in the BINDER quality standards.

11.3 Sending the chamber back to the manufacturer / distributor / service center

If you return a product to us for repair or any other reason, we will only accept the product upon presentation of an **authorization number** (RMA number) that has previously been issued to you. An authorization number will be issued after receiving your complaint either in writing or by telephone **prior** to your sending the BINDER product back to us. The authorization number will be issued following receipt of the information below:

- Product type and serial number
- Date of purchase
- Name and address of the dealer from which you bought the BINDER product
- Exact description of the defect or fault
- Complete address, contact person and availability of that person
- Exact location of the BINDER product in your facility
- A contamination clearance certificate (chap. 16) must be faxed / emailed in advance

The authorization number must be applied to the packaging in such a way that it can be easily recognized or be recorded clearly in the delivery documents.

For safety reasons we cannot accept a chamber delivery if it does not carry an authorization number.

12. End-of-life disposal

12.1 Disposal of the transport packing

Packing element	Material	Disposal
Straps to fix packing on pallet	Plastic	Plastic recycling
Wooden transport box (option) with metal screws	Non-wood (compressed wood, IPPC standard)	Wood recycling
	Metal	Metal recycling
Pallet (from size 115 on)	Compressed wood (IPPC standard)	Wood recycling
Transport box with metal clamps	Cardboard	Paper recycling
	Metal	Metal recycling
Top cover (size 720)	Cardboard	Paper recycling
Edge protection	Styropor® or PE foam	Plastic recycling
Protection of doors and hinges	PE foam	Plastic recycling
Bag for operating manual	PE foil	Plastic recycling
Insulating air cushion foil (packing of optional accessories)	PE foil	Plastic recycling

If recycling is not possible, all packing parts can also be disposed of with normal waste.

12.2 Decommissioning

- ED720, FED720: Turn off the chambers at the main power switch (chap. 2.3).
- Disconnect the chamber from the power supply (pull the power plug).
- With option inert gas connection: Close the inert gas supply and remove the gas connection.

	WARNING
	High concentration of inert gas. Death by suffocation. ➤ Observe the relevant regulations for handling inert gases. ➤ Close the gas supply when decommissioning the chamber.

Temporal decommissioning: See indications for appropriate storage, chap.3.3.

- Final decommissioning: Dispose of the chamber as described in chap.12.2.

12.3 Disposal of the device on the territory of the Russian Federation

	NOTICE Impact on the environment ➤ For final decommissioning and disposal of the oven, please contact BINDER Service. ➤ Study the disposal rules described in the legislation on environmental protection.
---	---

Upon reaching the service life limit or due to failure, sterilization medical equipment must be decommissioned and disposed of in accordance with applicable law.

Waste hazard class: A (Epidemiologically safe, non-toxic waste, which is close to solid domestic waste in composition, not in contact with biological fluids or infectious patients).

Disposal of medical devices is carried out in accordance with Sanitary Rules and Regulations (SanPiN) 2.1.3684-21.

Before handing over the equipment for disposal, the equipment must be disinfected, both the sterilization chamber and the external parts.

Equipment should not be disposed of in household waste collection points.

Disposal and disposal works may only be carried out by organizations licensed to dispose of medical equipment.

When replacing a used high-efficiency filter, it is subject to mandatory subsequent.

13. Technical description


13.1 Factory calibration and adjustment

The ovens were calibrated and adjusted in the factory. All test equipment used is subject to the administration of measurement and test equipment. They are controlled and calibrated to a DKD-Standard at regular intervals. Calibration and adjustment were performed using standardized test instructions, according to the QM DIN EN ISO 9001 system applied by BINDER (certified since December 1996 by TÜV CERT).

13.2 Definition of usable volume

Useful volume is the space in which the technical characteristics of the oven are guaranteed.

The usable volume illustrated below is calculated as follows:

A, B, C = Internal dimensions (W, H, D)

a, b, c = Wall clearances

$$a = 0.1 \times A$$

$$b = 0.1 \times B$$

$$c = 0.1 \times C$$

$$V_{USE} = (A - 2a) \times (B - 2b) \times (C - 2c)$$

Figure: Determination of the useable volume

To determine the useful volume, be guided by the data from the technical characteristics of the medical device (Chapter 13.4-13.6).

But if the temperature change is not important to the consumer, he can load the entire internal chamber of the oven.

Do NOT place charging material outside the usable volume.

Do NOT load the usable volume by more than half to enable sufficient airflow inside the chamber.

Do NOT divide the usable volume into separate parts with large area samples.

Do NOT place samples too close to each other in order to permit circulation between them and thus obtain a homogenous distribution of temperature.

13.3 Over current protection

Single-phase devices are protected by one or two miniature fuses against over current, accessible from the outside. The miniature fuses are located at the rear of the chamber above the power cable connection. Each fuse holder is equipped with a fuse clip 5mm x 20 mm. A fuse may be replaced only with a substitute of the same ratings. Refer to the technical data of the respective device type.

Two-phase devices are equipped with a resettable miniature circuit breaker (combination element).

Three-phase devices are equipped with an internal miniature circuit breaker.

13.4 Technical data for natural convection ovens

Model		E28	ED23	ED53	ED115	ED240	ED400	ED720
Case dimensions								
Width, net	mm	580	435	560	710	815	1235	1165
Height, gross (including legs)	mm	405	495	625	735	965	1025	1590
Depth, net	mm	425	520	565	605	760	765	816
Depth, gross (including door handle and ventilation pipe)	mm	477	625	640	680	815	870	870
Ventilation pipe, outer diameter	mm	52	52	52	52	52	52	52
Doors								
Number of doors	pcs	1	1	1	1	1	2	2
Chamber dimensions								
Width	mm	400	222	360	510	610	1000	960
Height	mm	280	330	420	530	760	800	1280
Depth	mm	250	277	380	425	550	500	610
Internal volume	l	28	20	57	114	255	400	749
Racks								
Number of racks (series)	pcs	2	2	2	2	2	2	2
Number of racks (max.)	pcs	4	4	4	5	8	10	16
Load per rack, max (chrome, stainless steel rack)	kg	10	12	30	30	40	35	45
Total permitted load	kg	25	25	70	150	270	90	315
Racks dimensions								
Chrome racks, WxDxH	mm	398x245x42	220x270x54	359x340x47	509x380x47	609x510x47	999x495x54	959x565x47
Stainless steel racks, WxDxH	mm	---	220x270x54	359x340x47	509x380x47	609x510x47	999x495x54	959x565x47

Perforated stainless steel racks, WxDxH		mm	399x242x16	220x245x15	359x339x16	509x379x16	609x509x16	999x489x16	959x569x16
Weight									
Oven (empty)		kg	22	26	39	54	85	125	161
Chrome racks		kg	0,628	0,483	0,794	1,536	0,352	6,037	5,043
Stainless steel racks		kg	-	0,483	0,799	1,546	2,367	6,021	5,081
Perforated stainless steel racks		kg	0,544	0,485	0,934	1,420	2,194	3,642	4,366
Temperature characteristics									
Temperature range	from	°C	60	---	---	---	---	---	---
	from ... degrees above ambient temperature	°C	---	5	5	5	5	5	5
	to	°C	230	300	300	300	300	300	300
Temperature fluctuation at 70°C		± K	1,5	0,2	---	---	---	0,2	---
Temperature fluctuation at 150°C		± K	---	0,5	0,4	0,4	1,0	0,5	0,8
Temperature variation at 70°C		± K	3	1,5	---	---	---	1,7	---
Temperature variation at 150°C		± K	---	2,8	2,5	2,0	2,0	3,0	3,2
Heating time up to 70°C		min	28	22	---	---	---	80	---
Heating time up to 150 °C		min	36	28	50	45	55	71	85
Heating time up to 250 °C		min	---	40	---	---	---	---	---
Heating time up to 300 °C		min	---	---	---	---	---	130	---
Temperature recovery time after opening the door for 30 seconds at 70°C		min	11	20	---	---	---	12	---
Temperature recovery time after opening the door for 30 seconds at 150°C		min	19	28	19	16	20	31	25
Temperature recovery time after opening the door for 30 seconds at 300°C		min	---	15	---	---	---	32	---
Ventilation data									
Ventilation rate at 70 °C		l/h	---	10	---	---	---	11	---
Ventilation rate at 150 °C		l/h	---	13	---	---	---	10	---
Ventilation rate at 300 °C		l/h	---	17	---	---	---	9	---
Electrical characteristics									
IP protection system according to 60529			IP20	IP20	IP20	IP20	IP20	IP20	IP20
Rated voltage (+/- 10%)	at a frequency of 50 Hz	V	230	230	230	230	230	400	400
	at a frequency of 60 Hz	V	230	230	230	230	230	400	400

Current type		1N~	1N~	1N~	1N~	1N~	3N~	3N~
Rated power	kW	0,8	0,8	1,05	1,25	2,25	3,4	4,10
Power cable length, (+/- 50 mm)	mm	1800	1800	2000	2000	2000	1800	1800
Power cable weight	kg	0,280	0,380	0,521	0,521	0,450	0,686	0,421
Power cord network connector		Grounding plug					CEE plug 5-pole	Grounding Plug
Fuse 5 x 20 mm / 250 V / delay or response time (T) external	A	6.3	10	6.3	6.3	12.5	---	---
Miniature circuit breaker (built-in)	A	---	---	---	---	---	16	16
Maximum temperature of independent protection system class 1	°C	320	---	330	330	320	---	320
Maximum temperature of independent protection system class 2	°C	---	320	300	300	300	320	300
Maximum temperature of independent protection system class 3.1	°C	---	350	300	300	300	350	300
Overvoltage category according to IEC 61010-1		II	II	II	II	II	II	II
Pollution degree according to IEC 61010-1		2	2	2	2	2	2	2
EMC acc. with EN 61326-1		Class A	Class A	Class A	Class A	Class A	Class A	Class A
Electric shock protection class		I	I	I	I	I	I	I
Characteristics in relation to the environment								
Energy consumption at 70°C	Wh/h	---	43	---	---	---	201	---
Energy consumption at 150°C	Wh/h	---	148	180	250	370	672	700
Energy consumption at 300°C	Wh/h	---	450	---	---	---	1000	---
Noise level (average value)	dB(A)	---	< 45	---	---	---	< 45	---
Model		E28	ED23	ED53	ED115	ED240	ED400	ED720

Note:

All technical data are based on empty ovens with standard equipment at ambient temperatures up to +22°C +/-3°C and voltage fluctuations +/-10%.

The technical data is determined according to the BINDER works standard Part 2:2015 and DIN 12880:2007.

The permissible deviation for the indicated dimensions in the finished product is +/-5 mm, the diameter of the ventilation pipe is +/-0.5 mm.

The calculated chamber volume is rounded to a whole number taking into account any specific features of the chamber internal surface.

The permissible deviation of the product weight is 5%.

All readings are average values and are typical for ovens manufactured in this series. We reserve the right to change technical specifications at any time.

If the oven chamber is fully loaded, the specific heating times may vary depending on the load.

13.5 Technical characteristics of ovens with forced convection.

Model		FD23	FD53	FD115	FD240
Case dimensions					
Width, net	mm	435	560	710	815
Height, gross (including legs)	mm	495	625	735	965
Depth, net	mm	520	565	605	760
Depth, gross (including door handle and ventilation pipe)	mm	616	640	680	815
Ventilation pipe, outer diameter	mm	52	52	52	52
Doors					
Number of doors	pcs	1	1	1	1
Chamber dimensions					
Width	mm	222	400	550	650
Height	mm	330	440	550	780
Depth	mm	277	345	385	512
Internal volume	l	20	60	116	259
Racks					
Number of racks (series)	pcs	2	2	2	2
Number of racks (max.)	pcs	4	4	5	8
Load per rack, max (chrome, stainless steel rack)	kg	12	30	30	40
Load per rack, max (perforated rack)	kg	--	35	35	35
Total permitted load	kg	25	70	150	270
Racks dimensions					
Chrome racks, WxDxH	mm	220x270x54	399 x 340 x 47	549 x 380 x 47	649 x 510 x 47
Stainless steel racks, WxDxH	mm	220x270x54	399 x 340 x 47	549 x 380 x 47	649 x 510 x 47
Perforated stainless steel racks, WxDxH	mm	220x245x15	399 x 319 x 16	549 x 359 x 16	649 x 489 x 16
Weight					
Oven (empty)	kg	27	39	54	85
Chrome racks	kg	0,483	0,851	1,633	2,481
Stainless steel racks racks	kg	0,483	0,856	1,643	2,497
Perforated stainless steel racks	kg	0,485	0,963	1,423	2,235
Temperature characteristics					
Temperature range	from ... degrees above ambient temperature	°C	5	10	10
	to	°C	300	300	300
Temperature fluctuation at 150°C	± K	0,3	0,3	0,3	0,5
Temperature variation at 150°C	± K	2,5	1.7	1.7	1.9
Heating time up to 150°C	min	25	15	19	20
Heating time up to 250°C	min	33	30	40	35
Temperature recovery time after opening the door for 30 seconds at 150°C	min	6	4	5	6
Ventilation data					
Ventilation rate	at 100°C	l/h	-	80	32
Ventilation rate	at 150°C	l/h	64	---	---
Electrical characteristics					
IP protection system according to EN 60529		IP20	IP20	IP20	IP20
Rated voltage (+/- 10%)	at a frequency of 50 Hz	V	230	230	230

	at a frequency of 60 Hz	V	230	230	230	230
Current type			1N~	1N~	1N~	1N~
Rated power	kW		0,8	1,10	1,30	2,30
Power cable length, (+/-50 mm)	mm		1800	2000	2000	2000
Power cable weight	kg		0,380	0,521	0,521	0,521
Power cord network connector			Grounding plug			
Fuse 5 x 20 mm / 250 V / delay or response time (T) external	A		10	6.3	6.3	12.5
Maximum temperature of independent protection system class 1	°C		---	330	350	350
Maximum temperature of independent protection system class 2	°C		320	300	300	300
Maximum temperature of independent protection system class 3.1	°C		350	300	300	300
Overvoltage category according to IEC 61010-1			II	II	II	II
Pollution degree according to IEC 61010-1			2	2	2	2
EMC acc. with EN 61326-1			Class A	Class A	Class A	Class A
Electric shock protection class			I	I	I	I
Characteristics in relation to the environment						
Noise level (average)	dB (A)		-	43	43	43
Energy consumption at 150°C	Wh/h		300	290	340	420
Model			FD23	FD53	FD115	FD240

Note:

All technical data are based on empty ovens with standard equipment at ambient temperatures up to +22°C +/-3°C and voltage fluctuations +/-10%.

The technical data is determined according to the BINDER works standard Part 2:2015 and DIN 12880:2007.

The permissible deviation for the indicated dimensions in the finished product is +/-5 mm, the diameter of the ventilation pipe is +/-0.5 mm.

The calculated chamber volume is rounded to a whole number taking into account any specific features of the chamber internal surface.

The permissible deviation of the product weight is 5%.

All readings are average values and are typical for ovens manufactured in this series. We reserve the right to change technical specifications at any time.

If the oven chamber is fully loaded, the specific heating times may vary depending on the load.

13.6 Technical data for natural and forced convection ovens

Model		FED53	FED115	FED240	FED400	FED720
Case dimensions						
Width, net	mm	560	710	815	1235	1165
Height, gross (including legs)	mm	625	735	965	1025	1590
Depth, net	mm	565	605	760	765	816
Depth, gross (including door handle and ventilation pipe)	mm	640	680	815	875	870
Ventilation pipe, outer diameter	mm	52	52	52	52	52
Doors						
Number of doors	pcs	1	1	1	2	2
Chamber dimensions						

Best conditions for your success							
Width	mm	400	550	650	1000	1000	
Height	mm	440	550	780	800	1300	
Depth	mm	345	385	512	500	570	
Internal volume	l	60	116	259	400	741	
Racks							
Number of racks (series)	pcs	2	2	2	2	2	
Number of racks (max.)	pcs	4	5	8	10	16	
Load per rack, max (chrome, stainless steel rack)	kg	30	30	40	35	45	
Load per rack, max (perforated rack)	kg	35	35	35	---	35	
Total permitted load	kg	70	150	270	90	315	
Racks dimensions							
Chrome racks, WxDxH	mm	399 x 340 x 47	549 x 380 x 47	649 x 510 x 47	999 x 495 x 54	999 x 570 x 47	
Stainless steel racks, WxDxH	mm	399 x 340 x 47	549 x 380 x 47	649 x 510 x 47	999 x 495 x 54	999 x 570 x 47	
Perforated stainless steel racks, WxDxH	mm	399 x 319 x 16	549 x 359 x 16	649 x 489 x 16	999 x 489 x 16	999 x 549 x 16	
Weight							
Oven (empty)	kg	39	54	85	145	162	
Chrome racks	kg	0,851	1,633	2,481	6,037	5,230	
Stainless steel racks racks	kg	0,856	1,643	2,497	6,021	5,269	
Perforated stainless steel racks	kg	0,963	1,423	2,235	3,642	4,345	
Temperature characteristics							
Temperature range	from ... degrees above ambient temperature	°C	10	10	10	5	12
	to	°C	300	300	300	300	300
Temperature fluctuation at 150°C	± K	0,3	0,3	0,5	0,7	0,6	
Temperature variation at 150°C	± K	1,4	1,2	1,6	3,8	2,0	
Heating time up to 150 °C	min	15	19	20	29	25	
Heating time up to 250 °C	min	30	40	35	55	40	
Temperature recovery time after opening the door for 30 seconds at 150 °C	min	4	5	6	6	6	
Ventilation data							
Ventilation rate	at 100 °C	l/h	80	32	9	---	5
Ventilation rate	at 150°C	l/h	---	---	---	18	---
Electrical characteristics							
IP protection system according to EN 60529		IP20	IP20	IP20	IP20	IP20	
Rated voltage (+/-10%)	at a frequency of 50 Hz	V	230	230	230	400	400
	at a frequency of 60 Hz	V	230	230	230	400	400
Current type			1N~	1N~	1N~	3N~	3N~
Rated power	kW	1,10	1,30	2,30	3,4	4,50	
Power cable length, (+/-50 mm)	mm	2000	2000	1800	1800	1800	
Power cable weight	kg	0,521	0,521	0,440	0,686	0,421	
Power cord network connector	Grounding plug						
Fuse 5 x 20 mm / 250 V / delay or response time (T) external	A	6.3	6.3	12.5	---	---	
Miniature circuit breaker (built-	A	---	---	---	16	16	

in)						
Maximum temperature of independent protection system class 1	°C	330	350	350	---	350
Maximum temperature of independent protection system class 2	°C	300	300	300	320	300
Maximum temperature of independent protection system class 3.1	°C	300	300	300	350	300
Overvoltage category according to IEC 61010-1		II	II	II	II	II
Pollution degree according to IEC 61010-1		2	2	2	2	2
EMC acc. with EN 61326-1		Class A	Class A	Class A	Class A	Class A
Electric shock protection class		I	I	I	I	I
Characteristics in relation to the environment						
Noise level (average)	dB (A)	43	43	43	-	43
Energy consumption at 150°C	Wh/h	300	340	420	1200	800
Model		FED53	FED115	FED240	FED400	FED720

Note:

All technical data are based on empty ovens with standard equipment at ambient temperatures up to +22°C +/-3°C and voltage fluctuations +/-10%.

The technical data is determined according to the BINDER works standard Part 2:2015 and DIN 12880:2007.

The permissible deviation for the indicated dimensions in the finished product is +/-5 mm, the diameter of the ventilation pipe is +/-0.5 mm.

The calculated chamber volume is rounded to a whole number taking into account any specific features of the chamber internal surface.

The permissible deviation of the product weight is 5%.

All readings are average values and are typical for ovens manufactured in this series. We reserve the right to change technical specifications at any time.

If the oven chamber is fully loaded, the specific heating times may vary depending on the load.

13.7 Overview of oven equipment and functions

To operate the oven, use only original BINDER accessories or accessories from third-party suppliers authorized by BINDER. The user is responsible for any risk caused by using unauthorized accessories.

The table below defines the composition of standard equipment and functions installed in a medical device by model.

√ - option availability

--- - lack of option

Description	E28
Hydraulic-mechanical thermostat	√
Control thermometer	√
Mechanical Timer 0-120 minutes	√
Independent thermal protection system Class 1	√

Adjustable ventilation flap

Description	ED23 FD23 ED400	FED400	ED53 ED115 ED240 ED720	FD53 FD115 FD240	FED53 FED115 FED240 FED720
Microprocessor temperature controller with timer	√	√	√	√	√
RS 422 interface	---	√	---	---	---
USB connector to read out the measured values	---	---	√	√	√
Ethernet interface	---	---	---	---	√
Timer functions: Delayed On, Delayed Off and Temperature dependent Delayed Off	---	√	---	---	√
Timer function: Delayed Off	√	---	√	√	---
Independent thermal protection system Class 1	---	---	√	√	√
Independent thermal protection system Class 2	√	√	√	√	√
Independent thermal protection system Class 3.1	√	√	√	√	√
Adjustable slow/linear heating function	√	√	√	√	√
Ventilation connection on the rear wall, with adjustable ventilation flap	√	√	√	√	√
Adjustable air exchange by means of a ventilation pipe	√	√	√	√	√

13.8 Accessories, spare parts and additional services (extract) - articles

	BINDER GmbH is only responsible for the safety of its devices if experienced electricians or qualified personnel authorized by BINDER have carried out all maintenance and repairs, and if components related to the safety of the oven are replaced with original spare parts in the event of a breakdown. The user is responsible for any risk caused by using unauthorized accessories.
---	--

13.8.1 Spare parts

Size	E28	ED23 FD23	ED53 FD53 FED53	ED115 FD115 FED115	ED240 FD240 FED240	ED400 FED400	ED720 FED720
Description	Art.No.						
Chromed rack for ED ovens	6004-0001	6004-0050	6004-0174	6004-0175	6004-0177	6004-0005	6004-0179
Chromed rack for FD, FED ovens	---	6004-0050	6004-0166	6004-0167	6004-0169	6004-0005	6004-0171
Stainless steel rack for ED ovens	---	6004-0051	6004-0158	6004-0159	6004-0161	6004-0011	6004-0163
Stainless steel rack for FD, FED ovens	---	6004-0051	6004-0150	6004-0151	6004-0153	6004-0011	6004-0155
Perforated stainless steel rack for ED ovens	6004-0028	6004-0052	6004-0190	6004-0191	6004-0193	6004-0032	6004-0195
Perforated stainless steel rack for FD, FED ovens	---	6004-0052	6004-0182	6004-0183	6004-0185	6004-0032	6004-0187
Replaceable HEPA fresh-air filters, class 14 (DIN EN 1822:2009)	---	---	6014-0003				

(for FD and FED only)							
Silicone door seal	6005-0001	6005-0090	6005-0254	6005-0255	6005-0258	6005-0069	6005-0260
Fluoroelastomer (FKM) door gasket (heat resistant up to 200°C)	---	8012-0493	6005-0265	6005-0266	6005-0268	8012-0497	6005-0270
Fuse 5 x 20 mm / 250 V / 6.3 A / delay or response time (T)	5006-0092	---	5006-0092	5006-0092	---	---	---
Fuse 5x20mm 250V 10A (M)	---	5006-0012	---	---	---	---	---
Fuse 5 x 20 mm / 250 V / 12.5 A / delay or response time (T)	---	---	---	---	5006-0096	---	---
Thermometer	5016-0002	---	---	---	---	---	---
Round timer switch 0-120 minutes	5017-0014	---	---	---	---	---	---

13.8.2 Additional services from BINDER.

For services not listed here, please contact BINDER Service.

Model	ED53 ED115 ED240 ED720	FD53 FD115 FD240	FED53 FED115 FED240 FED720
Installing the Ethernet communication Interface	√	√	---


Model	ED23	FD23	ED400
Installing the RS 422 communication Interface	√	---	√

Model			ED	FD	FED
Validation procedures	Art.No.				
IQ-OQ qualification folder			8012-0761	8012-0792	8012-0855
IQ-OQ-PQ qualification folder			8012-0942	8012-0943	8012-0944
Conducting IQ-OQ			DL400100	DL400100	DL400100
Conducting IQ-OQ-PQ			DL440500	DL440500	DL440500


Model			ED	FD	FED
Calibration services	Art.No.				
Temperature measurement calibration certificate (1 point)			DL300101	DL300101	DL300101
Certificate of calibration of spatial temperature measurement (at 9 points)			DL300109	DL300109	DL300109
Certificate of calibration of spatial temperature measurement (at 27 points)			DL300127	DL300127	DL300127
Measurement of air exchange rate according to ASTM D 5374, including certificate			DL330000	DL330000	DL330000

13.9 Medical device dimensions


13.9.1 Dimensions of oven E28 (mm)



13.9.2 Dimensions of ovens ED23, FD23 (mm)

13.9.3 Dimensions of ovens ED53, FD53, FED53 (mm)






13.9.4 Dimensions of ovens ED115, FD115, FED115 (mm)


ED



FD, FED

13.9.5 Dimensions of ovens ED240, FD240, FED240 (mm)



13.9.6 Dimensions of ovens ED400, FED400 (mm)


ED

FED

13.9.7 Dimensions of ovens ED720, FED720 (mm)

14. Certificates and declarations of conformity

CE	EU-Konformitätserklärung / EU Declaration of Conformity / Déclaration de conformité UE / Declaración de conformidad UE / Dichiarazione di conformità UE / Декларация соответствия EU
-----------	--

Hersteller / Manufacturer / Fabricant / Fabricante / Fabbricante / Производитель	BINDER GmbH БИНДЕР ГмбХ
Anschrift / Address / Adresse / Dirección / Indirizzo / Адрес	Im Mittleren Ösch 5, 78532 Tuttlingen, Germany Им Миттлерен Ош 5, 78532 Туттлинген, Германия

Produkt / Product / Produit / Producto / Prodotto / Продукт	Trocken- und Wärmeschränke mit freier Konvektion Drying and heating ovens with natural convection Etuves de chauffage et de séchage à convection naturelle Estufas de secado y calentamiento de convección natural Stufe per essiccazione e riscaldamento a convezione naturale Суховоздушные шкафы-стерилизаторы с естественной конвекцией
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	ED 23; ED 400

Produkt / Product / Produit / Producto / Prodotto / Продукт	Trocken- und Wärmeschränke mit Umluft Drying and heating ovens with forced convection Etuves de chauffage et de séchage à convection forcée Estufas de secado y calentamiento de convección forzada Stufe per essiccazione e riscaldamento a convezione forzata Суховоздушные шкафы-стерилизаторы с принудительной конвекцией
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	FD 23

Produkt / Product / Produit / Producto / Prodotto / Продукт	Trocken- und Wärmeschränke mit freier Konvektion und Umluft Drying and heating ovens with natural and forced convection Etuves de chauffage et de séchage à convection naturelle et forcée Estufas de secado y calentamiento de convección natural y forzada Stufe per essiccazione e riscaldamento a convezione naturale e forzata Суховоздушные шкафы-стерилизаторы с естественной и принудительной конвекцией
Typenbezeichnung / Type / Type / Tipo / Tipo / Тип	FED 400

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info@binder-world.com
www.binder-world.com
UK Authorized Representative:
Comply Express Ltd, Unit C2,
Coalport House, Stafford Park 1,
Telford TF3 3BD

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The machines described above are in conformity with the following EC/EU Directives (as published in the Official Journal of the European Union):

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The machines described above are conform to the mentioned EC/EU directives in regard to the relevant safety and health demands due to their conception and style of construction as well as to the version put onto market by us.

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The machines described above, corresponding to this, bear the CE-mark.

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

info@binder-world.com
 www.binder-world.com
 UK Authorised Representative:
 Comply Express Ltd, Unit C2,
 Coalport House, Stafford Park 1,
 Telford TF3 3BD

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info@binder-world.com
www.binder-world.com
UK Authorized Representative:
Comply Express Ltd, Unit C2,
Cowport House, Stafford Park 1,
Telford TF3 3BD

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Handwritten signature of J. Bollaender

J. Bollaender
Leiter F & E
Director R & D
Chef de service R&D
Responsable I & D
Direttore R & D
Глава департамента R&D

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info@binder-world.com
www.binder-world.com
UK Authorised Representative:
Comply Express Ltd, Unit C1,
Coalport House, Stafford Park 1,
Telford TF3 3BB

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0024177

15. Product registration

Online registration of equipment is possible at the link: www.binder-world.com/register

Easy registration in 3 steps:

1. Find the serial number of the product on the nameplate
2. Go to www.binder-world.com/register
3. Register your product by serial number

Online Product Registration

Register your BINDER now!

www.binder-world.com/register

The registration is free and takes just a few seconds

Advantages:

- ▶ Short response times if service is needed
- ▶ Fair prices when relocating or installing equipment
- ▶ Calibration as required at no charge in case of recalls
- ▶ Free information on news, product upgrades and accessories

Easy registered in 3 steps:

1. List serial number here: -
2. Go online: www.binder-world.com/register
3. Register serial number

16. Environmental cleanliness certificate

Declaration with regard to safety and health

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The German Ordinance on Hazardous Substances (GefStoffV), and the regulations regarding safety at the workplace, require that this form be filled out for all products that are returned to the manufacturer/distributor/service center, so that the safety and health of our employees can be warranted.

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ACT

On decontamination and decommissioning of equipment

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17. Manufacturer's warranties

The manufacturer, BINDER GmbH, guarantees that the device complies with the requirements of the technical documentation provided that the operating, storage, transportation and installation conditions established by the specified technical conditions and these operating instructions are observed.

The standard service life of the ovens is 10 years from the date of shipment by the manufacturer. The specified standards are provided on condition that the ovens are operated in accordance with the operating manual, the use of original spare parts for repairs and regular maintenance by service center specialists. The service life can be increased subject to regular maintenance based on the conclusion of a specialized organization.

Warranty period for dry-air ovens is 18 months from the date of putting the device into operation, but not more than 24 months from the date of shipment of the device by the manufacturer. The warranty period can be extended under a separate agreement.

Warranty repairs of ovens are carried out by supplier companies that have service departments authorized by the manufacturer. An approximate list of service centers authorized to carry out maintenance and repair of equipment is listed in Appendix 1, can be changed and expanded.

When carrying out warranty repairs, the warranty period is extended by the time elapsed from the moment the oven was received for repair until the end of the repair.

The consumer is deprived of the right to warranty repairs in the following cases:

- in case of damage to the oven during transportation in case of damage to the original packaging or its absence;
- in the presence of external damage after unpacking new equipment: case, controller unit, network cable and other components;
- the serial number tag or label has been removed from the device or component;
- if the factory seals or warranty stickers are damaged as a result of the consumer opening the case or equipment components;
- in case of damage caused by improper connection of equipment to utilities;
- in case of damage associated with violation of the operating rules given in the technical documentation for Drying and heating ovens.
- when the oven is out of order as a result of failures in the electrical network.
- in case of any malfunctions caused by external damage to the equipment, or excessive physical exertion during its transportation when the tilt or shock sensors are triggered.
- in the presence of mechanical damage to the external or internal parts, assemblies, conductors of the device, which arose during operation;
- if the device is damaged as a result of unqualified maintenance or repair: improper installation or damage to electronic components or mechanical parts.
- in case of any malfunctions caused by the operation of the equipment outside the declared range of controlled characteristics.
- in case of violation of the rules of transportation, storage and operation of the device;
- in case of damage caused by ingress of foreign substances, objects, liquids, insects and animals inside the device;
- in other cases, when the manufacturing enterprise has established that the consumer's claim is unreasonable.

In cases of oven breakdown in the post-warranty period, repairs can be carried out by service centers under a separate agreement at the expense of the consumer.

When contacting the service center, be prepared to provide the following information:

- serial number of the device;
- the date of receipt of the device from the supplier or dealer and the number of the document by which it was received;
- date of commissioning;
- a description of the external manifestation of the breakdown;
- what components and parts are broken, worn out, etc.

Annex 1.

LIST OF SERVICE CENTERS FOR REPAIR OF BINDER GMBH EQUIPMENT IN THE TERRITORY OF THE RUSSIAN FEDERATION.

ООО "DMA Service", 142180, Moskovskaya oblast, city of Podolsk, ulitsa Suvorova (Klimovsk microdistrict), Building 53a, office 2, tel. (495) 369 51 24; dma@dmaserv.ru

ООО "Vita-Pool", 125212, Moscow, Vyborgskaya ulitsa, house 16, building 1, office XVI, room 36, tel. (495) 514 19 07; service@vitapool.ru

ООО "Millab", 127410, Moscow, Inzhenernaya ulitsa, 18-1, 43, tel. (499) 933 71 47; service@millab.ru

ООО "Diaem", 129345, Moscow, Magadanskaya ulitsa, house 7, building 3, floor 2, office 210, tel. (495) 745 05 07; service@dia-m.ru

ООО NPP Labtech, 105043, Moscow, 3rd Parkovaya ulitsa, 8/19, room 9, 10, tel. (495) 276 77 00; service@labteh.ru

ООО "Helikon Company", 119619, Novomeshchersky proyezd, house 9, building 1, room 11, tel. (499) 705 50 50; service@helicon.ru

ООО "Amikton", 117452, Moscow, Chernomorsky Boulevard, house 17 building 1, floor/room/office 5/II/1/19, tel. (495) 236 72 64; service@amikton.ru

ООО "Nikolab", 115193, Moscow, ulitsa Petra Romanova, house 7, building 1, office I, room 16a, tel. (495) 675-26-67, service@nikolab.ru

For up-to-date information on service centers, you can contact the manufacturer, BINDER GmbH, <http://www.binder-world.com>, by e-mail east@binder-world.com, service@binder-world.com.

18. Information about the authorized representative on the territory of the Russian Federation

An authorized representative representing the interests of BINDER GmbH in the circulation of a medical device **Drying and heating ovens Binder in versions** on the territory of the Russian Federation, including on the issues of conformity assessment procedures and state registration, is:

ООО "European Medical Commonwealth", 117405, Moscow, ul. Dorozhnaya, 60B, www.medreg.ru, e-mail: med.emc@yandex.ru.

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The machines described above are in conformity with the following harmonized standards:

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
info@binder-world.com
www.binder-world.com
UK Authorised Representative:
Comply Express Ltd, Unit C2,
Caeport House, Stafford Park 1,
Telford TF3 3BD

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info@binder-world.com
www.binder-world.com
UK Authorised Representative:
Comply Express Ltd, Unit C2,
Coulport House, Stafford Park 1,
Telford TF3 3BD

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