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## CERTIFICATE

06446

中国国际贸易促进委员会暨中国国际商会  
China Council for the Promotion of International Trade is China Chamber of International Commerce

01807416

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China Council for the Promotion of International Trade  
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of Siemens Shanghai Medical Equipment Ltd. on the annexed DOCUMENT is genuine.

China Council for the Promotion  
of International Trade

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Authorized  
Signature:

Zhou Pingfan

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(Date: Jun. 06, 2023)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Name  
Department

Cheng Xiao  
SHS DI XP QT MEQ

**To the Federal Service for Surveillance  
in Healthcare of Russian Federation**

Mobile  
E-mail  
Date

+86 (21) 20606826  
cheng.xiao@siemens-healthineers.com  
Jun 02, 2023

### **Self-Truthfulness Declaration for Product Registration**

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**Product:** MULTIX Impact C

Herewith we, Siemens Shanghai Medical Equipment Ltd., in our responsibility as legal manufacturer would like to confirm that the following documents is true and correct copy of the original document:

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
# MULTIX Impact C

Operator Manual  
VA20

[siemens-healthineers.com](http://siemens-healthineers.com)

ООО ЦРПМ  [INFO@CIRMI.RU](mailto:INFO@CIRMI.RU)  [WWW.GOSZDRAVNADZOR.RU](http://WWW.GOSZDRAVNADZOR.RU)

**SIEMENS**  
Healthineers 

# Legend

|                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|---------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|            | Indicates a hint<br>Is used to provide information on how to avoid operating errors or information emphasizing important details                                                                                                                                                                                                                                                                                                              |
|            | Indicates the solution of a problem<br>Is used to provide troubleshooting information or answers to frequently asked questions                                                                                                                                                                                                                                                                                                                |
|            | Indicates a list item                                                                                                                                                                                                                                                                                                                                                                                                                         |
|            | Indicates a prerequisite<br>Is used for a condition that has to be fulfilled before starting a particular operation                                                                                                                                                                                                                                                                                                                           |
|            | Indicates a one-step operation                                                                                                                                                                                                                                                                                                                                                                                                                |
|            | Indicates steps within operating sequences                                                                                                                                                                                                                                                                                                                                                                                                    |
| <i>Italic</i>                                                                               | Is used for references and for table or figure titles                                                                                                                                                                                                                                                                                                                                                                                         |
|          | Is used to identify a link to related information as well as previous or next steps                                                                                                                                                                                                                                                                                                                                                           |
| <b>Bold</b>                                                                                 | Is used for screen outputs of the system                                                                                                                                                                                                                                                                                                                                                                                                      |
| Menu > Menu Item                                                                            | Is used for the navigation to a certain submenu entry                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Text</b>                                                                                 | Is used to identify inputs you need to provide                                                                                                                                                                                                                                                                                                                                                                                                |
| Orange                                                                                      | Is used to emphasize particularly important sections of the text                                                                                                                                                                                                                                                                                                                                                                              |
|  CAUTION | CAUTION<br>Used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury or material damage.<br>CAUTION consists of the following elements: <ul style="list-style-type: none"> <li>• Information about the nature of a hazardous situation</li> <li>• Consequences of not avoiding a hazardous situation</li> <li>• Methods of avoiding a hazardous situation</li> </ul> |
|  WARNING | WARNING<br>Indicates a hazardous situation which, if not avoided, could result in death or serious injury.<br>WARNING consists of the following elements: <ul style="list-style-type: none"> <li>• Information about the nature of a hazardous situation</li> <li>• Consequences of not avoiding a hazardous situation</li> <li>• Methods of avoiding a hazardous situation</li> </ul>                                                        |

|          |                                            |           |
|----------|--------------------------------------------|-----------|
| <b>1</b> | <b>Introduction</b>                        | <b>8</b>  |
| 1.1      | Intended purpose                           | 8         |
| 1.1.1    | Indications for use / Intended use         | 8         |
| 1.1.2    | Contraindications                          | 8         |
| 1.1.3    | Clinical benefits                          | 8         |
| 1.1.4    | Side effects                               | 9         |
| 1.1.5    | Residual risks                             | 9         |
| 1.1.6    | Patient group                              | 9         |
| 1.1.7    | Physical functionality                     | 9         |
| 1.1.8    | Intended users                             | 9         |
| 1.1.9    | Conditions of use                          | 9         |
| 1.1.10   | Minimum requirements concerning hardware   | 9         |
| 1.1.11   | IT network characteristics                 | 10        |
| 1.1.12   | Maintenance, cleaning, and disinfection    | 10        |
| 1.1.13   | Service                                    | 10        |
| 1.1.14   | Technical features                         | 10        |
| 1.1.15   | Frequently used operating functions        | 10        |
| 1.1.16   | Operating functions regarding safety       | 10        |
| 1.2      | Usability                                  | 10        |
| 1.3      | Information about this Operator Manual     | 11        |
| 1.3.1    | Area of application                        | 11        |
| 1.3.2    | Required documents                         | 11        |
| 1.3.3    | Operator Manuals on the internet           | 12        |
| 1.3.4    | Pictograms                                 | 12        |
| 1.3.5    | Names and parameters                       | 12        |
| 1.3.6    | Gender inclusivity                         | 13        |
| 1.3.7    | Illustrations                              | 13        |
| 1.3.8    | Value statements                           | 13        |
| <b>2</b> | <b>System Safety</b>                       | <b>14</b> |
| 2.1      | General safety information                 | 14        |
| 2.1.1    | User training                              | 14        |
| 2.1.2    | Laws and regulations                       | 14        |
| 2.1.3    | Explosion protection                       | 16        |
| 2.1.4    | Fire safety                                | 16        |
| 2.1.5    | Information about the software             | 17        |
| 2.1.6    | IT-security statement                      | 20        |
| 2.1.7    | Checks, tests                              | 20        |
| 2.1.8    | Maintenance plan                           | 21        |
| 2.1.9    | Installation, repair or modifications      | 26        |
| 2.1.10   | Error messages and status indicators       | 27        |
| 2.2      | Personal safety measures                   | 27        |
| 2.2.1    | Information about unit movements           | 27        |
| 2.2.2    | Danger zones / danger points               | 33        |
| 2.2.3    | Safety during patient examination          | 37        |
| 2.2.4    | Radiation protection                       | 40        |
| 2.3      | Equipment safety measures                  | 41        |
| 2.3.1    | Protection against electric shock          | 41        |
| 2.3.2    | Electromagnetic compatibility (EMC)        | 42        |
| 2.3.3    | Safety-relevant parts subject to wear      | 43        |
| 2.3.4    | Product service life                       | 43        |
| 2.3.5    | Mechanical safety                          | 43        |
| 2.3.6    | Combination with other products/components | 46        |
| 2.3.7    | Disposal                                   | 47        |

|          |                                                                                     |           |
|----------|-------------------------------------------------------------------------------------|-----------|
| 2.3.8    | Cleaning and disinfection                                                           | 47        |
| <b>3</b> | <b>System Overview</b>                                                              | <b>53</b> |
| 3.1      | Product description                                                                 | 53        |
| 3.1.1    | Application capabilities                                                            | 53        |
| 3.1.2    | System configurations                                                               | 54        |
| 3.2      | Basic system                                                                        | 54        |
| 3.2.1    | System overview in the examination room                                             | 54        |
| 3.2.2    | System overview in the control room                                                 | 55        |
| 3.2.3    | Image Station                                                                       | 55        |
| 3.2.4    | Generator ON/OFF console                                                            | 56        |
| 3.2.5    | Unit overview                                                                       | 57        |
| 3.3      | Ceiling-mounted tube support                                                        | 58        |
| 3.3.1    | Touch User Interface (TUI)                                                          | 59        |
| 3.3.2    | Multileaf collimator                                                                | 67        |
| 3.4      | Patient table                                                                       | 72        |
| 3.4.1    | Application                                                                         | 72        |
| 3.4.2    | Overview                                                                            | 72        |
| 3.4.3    | Brake knob                                                                          | 72        |
| 3.4.4    | Detector tray with wireless detector                                                | 73        |
| 3.4.5    | Emergency STOP push button                                                          | 74        |
| 3.4.6    | Foot kick switches in the table base                                                | 74        |
| 3.5      | Wall stand                                                                          | 75        |
| 3.5.1    | Description                                                                         | 75        |
| 3.5.2    | Application                                                                         | 75        |
| 3.5.3    | Operating elements                                                                  | 75        |
| 3.5.4    | Front panel markings                                                                | 77        |
| 3.6      | Wireless detectors (MAX wi-D or Core XL)                                            | 77        |
| 3.6.1    | Description                                                                         | 77        |
| 3.6.2    | Status and system messages                                                          | 79        |
| <b>4</b> | <b>System Operation</b>                                                             | <b>81</b> |
| 4.1      | On-Off                                                                              | 81        |
| 4.1.1    | Before using the system                                                             | 81        |
| 4.1.2    | Switching on a switched-off system                                                  | 81        |
| 4.1.3    | Switching on after power failure or after an emergency SHUTDOWN switch is activated | 82        |
| 4.1.4    | System switch-off                                                                   | 82        |
| 4.1.5    | Emergency SHUTDOWN switch activation                                                | 82        |
| 4.1.6    | Rebooting the system                                                                | 83        |
| 4.1.7    | UPS (for the imaging system)                                                        | 83        |
| 4.2      | Functional and safety check                                                         | 88        |
| 4.2.1    | Daily checks                                                                        | 88        |
| 4.2.2    | Monthly checks                                                                      | 91        |
| 4.2.3    | Functional check of the automatic systems                                           | 91        |
| 4.2.4    | Checks concerning radiology                                                         | 92        |
| 4.2.5    | DAP display test                                                                    | 92        |
| 4.3      | Preparation and movements of the overhead support                                   | 93        |
| 4.3.1    | Directions for moving and rotating the tube                                         | 93        |
| 4.3.2    | Directions for moving                                                               | 93        |
| 4.3.3    | Directions for rotating                                                             | 93        |
| 4.3.4    | Rotating the collimator                                                             | 93        |
| 4.3.5    | Positioning the 3D overhead support with tube unit                                  | 94        |
| 4.3.6    | Tube tracking                                                                       | 97        |

|          |                                                    |            |
|----------|----------------------------------------------------|------------|
| 4.3.7    | SID tracking                                       | 98         |
| 4.3.8    | Deactivation of a tracking control                 | 99         |
| 4.3.9    | SID display                                        | 99         |
| 4.3.10   | Rotating the multileaf collimator                  | 99         |
| 4.4      | Preparation and movements of the patient table     | 101        |
| 4.4.1    | Collision protection                               | 101        |
| 4.4.2    | Adjusting table height                             | 101        |
| 4.4.3    | Positioning the tabletop                           | 102        |
| 4.4.4    | Inserting wireless detector                        | 103        |
| 4.4.5    | Positioning the detector tray                      | 104        |
| 4.4.6    | Scattered radiation grid                           | 105        |
| 4.5      | Preparation and movements of the wall stand        | 108        |
| 4.5.1    | Safety information                                 | 108        |
| 4.5.2    | Detector unit movements                            | 109        |
| 4.5.3    | Inserting the detector into the wall stand         | 110        |
| 4.5.4    | Removing the detector                              | 111        |
| 4.5.5    | Scattered radiation grid                           | 112        |
| 4.5.6    | Positioning the tube unit automatically            | 114        |
| 4.5.7    | Positioning the tube unit manually                 | 114        |
| 4.5.8    | Centering the tube unit in vertical direction      | 115        |
| 4.5.9    | Moving the tube unit to the parking position       | 115        |
| 4.6      | Preparation and handling of the wireless detectors | 115        |
| 4.6.1    | General safety information                         | 115        |
| 4.6.2    | Recommendations                                    | 116        |
| 4.6.3    | Notes on handling the detector                     | 117        |
| 4.6.4    | Maintenance                                        | 120        |
| 4.6.5    | Operation and setting of WLAN                      | 120        |
| 4.6.6    | Battery change (MAX wi-D)                          | 122        |
| 4.6.7    | Battery change (Core XL)                           | 124        |
| 4.6.8    | Charging the battery of the detector               | 125        |
| 4.6.9    | FD management                                      | 126        |
| 4.6.10   | Attaching the detector                             | 127        |
| 4.6.11   | Detaching the detector                             | 128        |
| 4.6.12   | Detector sharing with MAX wi-D                     | 128        |
| 4.6.13   | Detector calibration                               | 128        |
| 4.7      | Emergency operation                                | 129        |
| <b>5</b> | <b>Examination</b>                                 | <b>130</b> |
| 5.1      | Safety information for examination                 | 130        |
| 5.2      | Positioning the patient on the table               | 131        |
| 5.2.1    | Hygiene                                            | 131        |
| 5.2.2    | Radiation protection                               | 131        |
| 5.2.3    | Grid and SID                                       | 131        |
| 5.2.4    | Table position                                     | 131        |
| 5.2.5    | Positioning very heavy patients                    | 131        |
| 5.2.6    | Tube unit position                                 | 132        |
| 5.2.7    | Manual SID                                         | 132        |
| 5.2.8    | Detector pulled-out                                | 132        |
| 5.3      | Positioning the patient at the wall stand          | 132        |
| 5.3.1    | Hygiene                                            | 132        |
| 5.3.2    | Radiation protection                               | 132        |
| 5.3.3    | Wireless detector                                  | 132        |
| 5.3.4    | Grid                                               | 133        |
| 5.3.5    | Overhead patient hand grip                         | 133        |

|          |                                                    |            |
|----------|----------------------------------------------------|------------|
| 5.3.6    | Tube unit position                                 | 133        |
| 5.3.7    | Operating Height                                   | 133        |
| 5.4      | IONTOMAT automatic exposure control                | 133        |
| 5.4.1    | Performing IONTOMAT automatic exposure control     | 134        |
| 5.5      | Additional Cu filter                               | 134        |
| 5.5.1    | Selecting additional filters                       | 135        |
| 5.6      | Light localizer and collimation                    | 135        |
| 5.6.1    | Switching on and off of the light localizer        | 136        |
| 5.7      | Collimation                                        | 137        |
| 5.7.1    | Collimating manually                               | 137        |
| 5.7.2    | Collimation / Automatic Collimation Size Sensing   | 137        |
| 5.7.3    | Restoring the last collimation                     | 138        |
| 5.8      | Exposures onto the patient table                   | 138        |
| 5.8.1    | Vertical projection onto the patient table         | 139        |
| 5.8.2    | Vertical projection onto the pulled-out detector   | 140        |
| 5.8.3    | Oblique projection onto the patient table          | 141        |
| 5.9      | Exposures onto the wall stand                      | 142        |
| 5.9.1    | Horizontal projection onto the wall stand          | 142        |
| 5.9.2    | Vertical projection onto the wall stand            | 143        |
| 5.9.3    | Oblique projection onto the wall stand             | 143        |
| 5.10     | Free exposure onto the wireless detector           | 144        |
| 5.10.1   | Preparations                                       | 145        |
| 5.10.2   | Positioning the tube unit                          | 145        |
| 5.10.3   | SID                                                | 145        |
| 5.10.4   | Positioning for vertical projection                | 146        |
| 5.10.5   | Positioning for oblique projection                 | 146        |
| 5.10.6   | Horizontal projection onto the mobile detector     | 146        |
| 5.11     | Pediatric Use: Summary                             | 147        |
| 5.11.1   | Introduction                                       | 147        |
| 5.11.2   | Use of anti-scatter grid                           | 147        |
| 5.11.3   | Use of filtration                                  | 147        |
| 5.11.4   | Automatic Exposure Control (AEC)                   | 148        |
| 5.11.5   | Specific organ programs for pediatric examinations | 148        |
| 5.11.6   | Pediatric positioning aids                         | 148        |
| 5.12     | Ortho operation                                    | 148        |
| 5.12.1   | Ortho on patient table                             | 149        |
| 5.12.2   | Ortho on wall stand                                | 149        |
| 5.12.3   | Organ program for Ortho                            | 150        |
| 5.12.4   | Preparing Ortho                                    | 150        |
| 5.12.5   | Releasing the Ortho exposures                      | 155        |
| <b>6</b> | <b>Accessories and Auxiliary Devices</b>           | <b>157</b> |
| 6.1      | Preliminary remarks                                | 157        |
| 6.1.1    | Equipment with accessories                         | 157        |
| 6.1.2    | Proper use of the product                          | 157        |
| 6.1.3    | Safety                                             | 158        |
| 6.1.4    | Orientation                                        | 159        |
| 6.2      | Description                                        | 159        |
| 6.2.1    | Trolley UM mobile patient table                    | 159        |
| 6.2.2    | Patient positioning mattress                       | 163        |
| 6.2.3    | Paper roll holder                                  | 164        |
| 6.2.4    | Hand grip strip                                    | 164        |
| 6.2.5    | Hand grip, lateral                                 | 166        |
| 6.2.6    | Overhead patient hand grip                         | 168        |

|          |                                                                                                        |            |
|----------|--------------------------------------------------------------------------------------------------------|------------|
| 6.2.7    | BABIX holder                                                                                           | 170        |
| 6.2.8    | Multipurpose / Ortho stand                                                                             | 171        |
| 6.2.9    | Compression Belt                                                                                       | 184        |
| 6.2.10   | Cassette holder, lateral                                                                               | 192        |
| 6.2.11   | Cassette holder without retainer                                                                       | 194        |
| 6.2.12   | Mobile detector holder CR                                                                              | 195        |
| 6.2.13   | Detector holder, lateral                                                                               | 196        |
| 6.2.14   | Clip-on grids                                                                                          | 198        |
| 6.2.15   | FD protector (patient support)                                                                         | 199        |
| 6.2.16   | Wall holder for grids                                                                                  | 200        |
| 6.2.17   | Compensation filters                                                                                   | 200        |
| 6.2.18   | Filter holder for eight filters                                                                        | 204        |
| 6.2.19   | Three-field template, set                                                                              | 205        |
| 6.2.20   | Wireless Remote Control Console                                                                        | 206        |
| 6.2.21   | Foot switch                                                                                            | 211        |
| 6.2.22   | Battery charger (MAX wl-D)                                                                             | 212        |
| 6.2.23   | Battery charger (Core XL)                                                                              | 214        |
| 6.2.24   | Wall charger                                                                                           | 215        |
| 6.2.25   | Barcode scanner                                                                                        | 216        |
| 6.2.26   | Intercom system in control room                                                                        | 217        |
| <b>7</b> | <b>Technical Description</b>                                                                           | <b>218</b> |
| 7.1      | Location of labels                                                                                     | 218        |
| 7.1.1    | Labels on the ceiling-mounted tube support                                                             | 218        |
| 7.1.2    | Labels on the tube unit with multileaf collimator                                                      | 221        |
| 7.1.3    | Labels on the wall stand                                                                               | 223        |
| 7.1.4    | Labels on the patient table                                                                            | 225        |
| 7.1.5    | Labels on the generator ON/OFF console                                                                 | 226        |
| 7.1.6    | Labels on the generator cover                                                                          | 227        |
| 7.2      | Technical data                                                                                         | 228        |
| 7.2.1    | System                                                                                                 | 228        |
| 7.2.2    | Components                                                                                             | 229        |
| 7.2.3    | Network information                                                                                    | 236        |
| 7.3      | Notice concerning electromagnetic compatibility (EMC)                                                  | 237        |
| 7.3.1    | Guidance and manufacturer's declaration - electromagnetic emissions                                    | 238        |
| 7.3.2    | Guidelines and manufacturer's declaration - electromagnetic immunity                                   | 238        |
| 7.3.3    | Recommended separation distance between portable and mobile RF communications equipment and the system | 240        |
| 7.3.4    | RF transmitters                                                                                        | 241        |
| <b>8</b> | <b>Glossary</b>                                                                                        | <b>242</b> |
|          | <b>Index</b>                                                                                           | <b>243</b> |

# 1 Introduction

To operate the MULTIX Impact C system accurately and safely, the operating personnel must have the necessary expertise as well as knowledge of the entire Operator Manual. It must be read carefully prior to starting up the MULTIX Impact C system.

## 1.1 Intended purpose

Intended Purpose Statement under the European Medical Device Regulation 2017/745:

Stationary X-ray system intended for various procedures that visualizes a variety of anatomical structures.

### 1.1.1 Indications for use / Intended use

MULTIX Impact C is a radiographic system used in hospitals, clinics, and medical practices. MULTIX Impact C enables radiographic exposures of the whole body including: skull, chest, abdomen, and extremities and may be used on pediatric, adult and bariatric patients. Exposures may be taken with the patient sitting, standing, or in the prone position. MULTIX Impact C is not intended for mammography.

MULTIX Impact C uses digital detectors for generating diagnostic images by converting X-rays into image signals. MULTIX Impact C is also designed to be used with conventional film/screen or Computed Radiography (CR) cassettes.

### 1.1.2 Contraindications

There are currently no known contraindications for MULTIX Impact C.

All residual risks are communicated to the user by warnings and cautions in the operator manual.

The final decision for use of the medical device in a certain application is made by the physician based on his/her medical knowledge and the risks of the case involved.

Clinically relevant hazards and risks were mitigated as summarized in the risk management report.

Clinically relevant side effects and risks that are generally observed with the diagnostic radiography and radiological technology and which are described in the literature, are the exposition of the patients to intended ionizing radiation.

### 1.1.3 Clinical benefits

The device is a digital X-ray system to generate X-ray images from the whole body including the skull, chest, abdomen and extremities. The acquired images support medical professionals to make diagnostic and/or therapeutic decisions. Generic clinical benefits of radiographic examinations within the intended use are applicable for this device.

### 1.1.4 Side effects

The side effects of the device are derived from the risk analysis. Based on the clinical evaluation, there are no known side effects associated with the use of MULTIX Impact C system.

### 1.1.5 Residual risks

Residual risks have been assessed and all relevant information is provided to the user via this Operator Manual.

Please note that risks inherently connected to the use of this device (e.g. radiation or moving parts) can not be avoided due to the nature of this device.

### 1.1.6 Patient group

The MULTIX Impact C system is designed for adult and pediatric (removable scattered radiation grid) patients.

### 1.1.7 Physical functionality

The MULTIX Impact C system may be used for projection of radiographic examinations of both adults and children (removable scattered radiation grid).

- On the patient table, X-rays can be taken on the recumbent or sitting patient from head to foot. Exposures in the region of the skull, of the spine, of thorax, lungs and abdomen as well as of the extremities are possible.
- On the wall stand, X-rays can be taken on the standing or sitting patient (chest, and so on).
- In addition, exposures onto the wireless detectors in the form of free exposures are possible.

### 1.1.8 Intended users

The MULTIX Impact C must be operated by qualified personnel with the necessary knowledge and expertise to operate X-ray equipment in accordance with the relevant country specific regulations.

(→ Page 14 *User training*)

### 1.1.9 Conditions of use

The system is suitable for medical facilities such as hospitals, doctor's offices, etc.

The climatic conditions specified in Technical data, Ambient conditions (operation) must be complied with.

### 1.1.10 Minimum requirements concerning hardware

The MULTIX Impact C system is being delivered, installed, and connected to the IT environment as a complete and functioning system by our service organization.

### 1.1.11 IT network characteristics

MULTIX Impact C may only be run in an environment approved or authorized by the manufacturer. The system may not be operated in public networks without firewalls or isolation.

Refer to the "MULTIX Impact C Security White Paper and MDS<sup>2</sup> Form".

### 1.1.12 Maintenance, cleaning, and disinfection

Refer to (→ Page 14 *System Safety*) and (→ Page 81 *System Operation*).

### 1.1.13 Service

The service must only be performed by qualified technical personnel.

### 1.1.14 Technical features

Refer to (→ Page 228 *Technical data*).

### 1.1.15 Frequently used operating functions

Refer to (→ Page 81 *System Operation*).

### 1.1.16 Operating functions regarding safety

Refer to (→ Page 14 *System Safety*).

## 1.2 Usability

The systems must only be used by persons with the necessary specialist knowledge after training, for example physicians, radiologists, cardiologists, and medical specialists.

- Patient population: newborn to geriatric.
- Operator profile: The usage of the system described in this Operator Manual requires specific technical and medical knowledge and skills regarding, at a minimum, radiation protection, safety procedures, and patient safety.  
  
People using, moving, working with the system must have acquired such knowledge and skills during their curriculum.
- Language understanding: Users must understand the language of the Operator Manual before touching the system.
- Equipment training: Application training is supplied with the equipment according to the hand-over contract. It is mandatory to follow such application training supplied by Siemens Healthineers Representative before any use of the system.

The follow-up training, which is necessary due to change of personnel, is in the responsibility of the operator of the system. Any additional training can be requested from Siemens Healthineers.

- **Operator Manual and precautions:** Read and understand all the instructions in the Operator Manual before trying to use the system and request additional training from Siemens Healthineers if necessary.

Keep the Operator Manual with the equipment at all times and periodically review the procedures and safety precautions.

Failure to follow the Operator Manual and safety precautions could result in serious injury to the patient, others or yourself.

- **Patient safety:** Give assistance for getting the patient on and off the table. Be sure all patient health lines (IV, oxygen, and so on) are positioned so they will not be caught when moving the equipment. Never leave the patient unattended while in the system room.

An unattended patient could fall from the table, start a movement control, or encounter other problems which could be hazardous.

- **Radiation safety:** Always use correct technique factors for each procedure to minimize X-ray exposure and to produce the best diagnostic results.
- **Establish emergency procedures:** It is not always possible to determine when some components, such as X-ray tubes, are nearing the end of their operating lives. These components could stop operating during a patient examination. Have an emergency plan in case of non-availability of the system, for example, to use a different system.

#### CAUTION

Operation of the system by untrained users

**Risk of incorrect diagnosis or treatment due to misinterpretation of image information**

- ◆ The system must only be used by persons with the necessary specialist knowledge, for example, physicians, trained radiologists, or trained technologists, after an appropriate application training.

## 1.3 Information about this Operator Manual

To improve readability, your complete Operator Manual has been broken down into several individual Operator Manuals with thematically distinct content:

- Operator Manual - MULTIX Impact C system
- Operator Manual - MULTIX Impact C imaging system

### 1.3.1 Area of application

This Operator Manual is valid for the following product:

- MULTIX Impact C (VA20)

### 1.3.2 Required documents

Reference is made in this Operator Manual to the following documents:

- Operator Manual - MULTIX Impact C imaging system

### 1.3.3 Operator Manuals on the internet

The electronic Operator Manuals for the system are available on the internet:  
<https://doclib.healthcare.siemens.com>.

### 1.3.4 Pictograms

The following are pictograms and their meanings as they may apply to your system (IEC standard).

|                                                                                     |                                            |
|-------------------------------------------------------------------------------------|--------------------------------------------|
|    | Alternating current                        |
|    | Caution: Read Operator Manual              |
|    | Caution: Laser                             |
|    | Pushing prohibited                         |
|   | Off (disconnection from line power supply) |
|  | On (connection to line voltage)            |
|  | Attention: Hot surface (at the X-ray tube) |
|  | Degree of protection type B                |

### 1.3.5 Names and parameters

All names and data on patients and equipment that are used as examples in this Operator Manual are entirely fictional.

Any resemblance to names of real persons and institutions is entirely coincidental.

All parameters and images shown in this Operator Manual are examples. Only the parameters displayed on your system are relevant.

The MULTIX Impact C may be shown with expanded equipment or with optional accessories in some figures.

### 1.3.6 Gender inclusivity

At Siemens Healthineers, we want to address all genders equally in our Instructions for Use. It is not our intention to exclude anyone. We continue to update our documentation with this in mind.

### 1.3.7 Illustrations

All illustrations of the equipment and of the user interface in this Operator Manual, are examples only.

Differences in detail may occur in your system due to the installed options, configuration and constant development of the system.

Reproduction of images can cause loss of detail. Pictures in this Operator Manual do not therefore provide any indication of image quality.

All names of patients in images or illustrations, are purely fictional. Any similarities with existing persons are entirely coincidental.

### 1.3.8 Value statements

All technical data are typical values unless specific tolerances are stated.

Values in pictures of the software user interface, have no clinical meaning.

Only set the values preset in the exam sets provided or the values recommended by experienced application specialists.

## 2 System Safety

### 2.1 General safety information

The organization responsible for the safety of the medical device according to IEC 60601 ff. is the user or the operator/operating company of the device.

#### 2.1.1 User training

To operate the MULTIX Impact C system accurately and safely, the operating personnel must have the necessary expertise as well as knowledge of the entire Operator Manual. It must be read carefully prior to starting up the MULTIX Impact C system.

Application training is supplied with the equipment according to the hand-over contract. It is mandatory to follow such application training supplied by Siemens Healthineers Representative before any use of the system.

The follow-up training, which is necessary due to change of personnel, is in the responsibility of the operator of the system. Any additional training can be requested from Siemens Healthineers.

#### CAUTION

Operation of the system by untrained users

**Risk of incorrect diagnosis or treatment due to misinterpretation of image information**

- ◆ The system must only be used by persons with the necessary specialist knowledge, for example, physicians, trained radiologists, or trained technologists, after an appropriate application training.

#### 2.1.2 Laws and regulations

The statutory regulations of the respective country must be followed.

If legally binding regulations govern the installation and/or the operation of the product, it is the responsibility of the installer or the operator to observe these regulations.

Modifying the medical product is not permitted.

#### Country-specific regulations

The legally established country-specific regulations must be followed in all countries. Deviating from this Operator Manual, values may be set according to country-specific regulations.

## Regulations in the EU

Due to the Medical Device Directive 93/42/EEC and EU regulation 2017/745 (MDR), MULTIX Impact C complies with the standards of the IEC 60601-1 series.

### Incident reporting

According to EU regulation 2017/745 (MDR), any serious incident that has occurred in relation to the device should be reported to the manufacturer and the competent authority of the EU Member State in which the user and/or patient is established.

## Additional national regulations

To ensure safety for operators, patients, and third parties, we recommend compliance with the regulations listed here to the extent that local national laws allow.

### Germany RöV (X-ray Ordinance) - Legally required tests according to §16:

#### Before start-up

- Acceptance test
- Expert inspection

#### During operation

- According to the time intervals specified

#### Constancy test

#### MedBetreibV - Ordinance for system owners of medical equipment

- Safety check for technical construction at least every two years

## Regulations in the USA and Canada

Corresponding to the Medical Device Directive in the EU, MULTIX Impact C complies in the U.S.A. and Canada with the applicable standards.

Federal law stipulates that this system may only be sold to a physician or by order of a physician.

This system is intended for use in radiographic examinations and under the guidance of a trained health care professional.

## Handling of non RoHS spare parts

In the interests of preserving natural resources, Siemens Healthineers may recycle spare parts which in some circumstances are not compliant with RoHS 2011/65/EU.

However, in these cases Siemens Healthineers uses the valid exemption to these regulations which allow the use of these parts.

## Data security

Personal data are subject to data protection.

- Please observe the relevant legal provisions.

### Legally required tests

Legally required tests must be performed at the specified intervals.

These tests include, for example,

- Constancy test according to the X-ray ordinance (§16 RöV) in the Federal Republic of Germany.
- Tests based on DHHS guidelines (Department of Health and Human Services) where applicable.

For more detailed information refer to (→ Page 88 *Functional and safety check*).

### 2.1.3 Explosion protection

MULTIX Impact C is not designed for operation in explosion-endangered areas.

#### CAUTION

The system is not designed for operation in explosion-endangered areas. It does not fulfill the requirements for AP/APG classification.

#### **Danger of explosion!**

- ◆ Do not use the system where explosive atmospheres can arise.

Only products with the AP mark (AP = anesthetics test) may be used in explosion-endangered areas

### 2.1.4 Fire safety

In the event of a fire, shut down the entire system immediately, that is, disconnect the system from the main power supply.

- 1 Press the emergency SHUTDOWN button or actuate the main or disconnecting switch.
- 2 Use a CO<sub>2</sub> fire extinguisher.

Do not use water!

**CAUTION**

Fire inside or in the vicinity of the system

**Risk of injury of patient and personnel and damage to property. Risk of inhalation of fumes from burning plastics.**

- ◆ Switch off the system in the event of fire.
- ◆ Help the patient to escape from the system.
- ◆ Make sure that you know where the escape routes are.
- ◆ Make sure that you know where the fire extinguishers are located and familiarize yourself with the use of them.

- 3 Contact Siemens Healthineers Customer Service before performing any refurbishing work and before starting the system up again.

## 2.1.5 Information about the software

### Network security

**CAUTION**

Unauthorized access to the system

**Hazards up to and including loss of patient data and non-operational system**

- ◆ The hospital is responsible for network security at the site.
- ◆ Set up user account password protections.
- ◆ Do not allow users to change configuration files.

Please note that the following are changes to the network that may introduce new risks:

- Network configuration change;
- Connection of additional systems or equipment to the network;
- Disconnection of systems or equipment from the network;
- Update or upgrade of equipment connected to the network.

Please note the following information relating to cybersecurity:

- Appropriate physical access controls shall be implemented;
- Cables shall be protected from interception, interference, and damage;
- Maintenance and inspection of devices shall only be carried out by authorized personnel;
- External personnel who perform maintenance or service work shall be accompanied at all times;
- Devices shall be secured with physical locking cables or securely locked away;
- The system shall lock down automatically or terminate the session after exceeding a reasonable defined idle time limit;

- Applications shall log off the user or terminate the session after exceeding a reasonable defined idle time limit;
- Malicious software shall not be created, downloaded, or distributed on purpose;
- User accounts shall only be assigned to and used by individual persons;
- For each account a different password shall be used;
- IP addresses assigned to systems shall be logged;
- Device log data shall be continuously analyzed;
- Device log data shall be stored for a reasonable time period;
- Devices shall be placed in network zones which fulfill and are appropriate for their information security requirements;
- The system must be used by persons with necessary specialist knowledge after training supplied by Siemens Healthineers Representative;
- The hospital is responsible for authorized access to the system;
- We recommend installing the system in a separate examination room with a door controller;
- The system can be connected to the internal network. The hospital is responsible for network security at the site;
- The Device Guard built into Windows 10 is used to prevent unauthenticated software from being executed in the system;
- The firewall built into Windows 10 is used to block network connections through illegal ports;
- The security service is used to encrypt/decrypt important files;
- The corespot service is used to enhance the security of the operating system;
- Normal users shall use a unique ID and password to log into the system;
- The system allows access in a special mode without authentication in order to enable examination in the event of an emergency. For this special mode no password is set;
- The system has basic requirements on the length/complexity/expiry time of the password;
- The system predefines several user roles, and users can only access the functions defined in their roles;
- The system has been tested before its launch, making sure that there are no software security vulnerabilities;
- The major 3rd party software used by the system has been registered to the Siemens Healthineers SVM service to monitor security vulnerabilities that may be found in future;
- Siemens Healthineers will analyze vulnerability notifications from SVM in a timely manner and take appropriate measures according to their severity;

- Siemens Healthineers will provide timely patches or mitigations for serious software security vulnerabilities;
- Siemens Healthineers releases a security patch for the system regularly (90 days) to repair general security vulnerabilities.

### Third-party software

Only software authorized by Siemens Healthineers for use with this product may be used.

The pre-installed security package provides protection against cyber attacks, viruses, malware and other damaging software. It ensures that only trusted applications are run on the systems. It blocks unauthorized access, provides protection from network threats and infected USB sticks and thus offers control over when and who may make changes.

#### CAUTION

Use of unapproved or manipulated software or hardware components

#### **Risk of malfunctions posing a hazard to patients and equipment**

- ◆ Only use software or hardware components approved by Siemens Healthineers.
- ◆ Perform repairs only with our expressed written approval.

### Original equipment manufacturer

Original Equipment Manufacturer (OEM), Open Source Software (OSS) from third-party providers in Siemens Healthineers products

#### Introduction

The product can contain OEM and OS software licensed and developed by third-party providers (depending on the system configuration and option).

These OEM and/or OS files are protected by copyright. Your right to use OEM and/or OS software beyond the mere execution of the Siemens Healthineers program is governed by the relevant terms of the OEM and/or OS software license.

Complying with these licensing conditions authorizes you to use the OEM and/or OS software as provided in the respective license. In case of conflicts between the licensing conditions of Siemens Healthineers and the OEM and/or OS software license conditions, the OEM and/or OS software conditions govern the OEM and/or OS portions of the software.

The OEM and/or OS software licenses (license conditions) are available on the CDs/DVDs (open source/legal concept) delivered together with the system.

#### For Open Source Software (OSS) only

Open source software is licensed free of charge (meaning that exercising the licensing rights is free of charge, while fees may be charged to recover the costs borne by Siemens Healthineers). If programs contained in this product are licensed under GNU General Public License (GPL), GNU Lesser General Public License (LGPL), or another open source license that requires source code, the corresponding code is delivered together with the device on CD/DVD.

Warranty regarding the further use of software from third-party providers:

Siemens Healthineers makes no warranty regarding the OEM and/or OS software programs in this device, when such programs are used in a manner other than the program execution intended by Siemens Healthineers. The license mentioned previously define the warranty, if any, from the authors or licensors of the OEM and/or OS software. Siemens Healthineers specifically disclaims any warranties for defects caused by altering any OEM and/or OS software program or the product's configuration. If the OEM and/or OS software violates the intellectual property of a third party, you have no warranty claim against Siemens Healthineers.

Technical support is available only for software that has not been modified.

## Copyright

The system and user software used in this product is copyright protected.

## DICOM conformity

The digital imaging system is DICOM-conform. The "DICOM Conformance Statement" is available on request from Siemens Healthcare or via the Siemens Healthineers Internet homepage.

## Remote software update

### CAUTION

Software installation will disturb system functions.

**Interruption of patient treatment and loss of patient data**

- ◆ Do not start installation during patient treatment.

### CAUTION

Installation of the Update package failed.

**Failed Software installation will result in an undefined system state.**

- ◆ Do not use the system and call the Siemens Healthineers Service.

## 2.1.6 IT-security statement

The device can be connected to an IT-network. In this case, detailed information regarding the combination and the risk caused by the combination are described in the IT-security statement of accompanying documents.

## 2.1.7 Checks, tests

Before the system is used for examinations, the user must ensure that all safety-relevant devices function correctly and that the system is ready for operation.

Important functional tests and checks must be made at certain time intervals.

(→ Page 88 *Daily checks*)

### 2.1.8 Maintenance plan

#### WARNING

Wear or material fatigue in the system or accessories

##### **Risk of injury or system damage**

- ◆ Follow the maintenance guidelines to maintain safety and proper function of the system.
- ◆ Check accessories for wear before use.

#### CAUTION

Image artifacts due to unperformed maintenance

##### **Incorrect diagnosis basis**

- ◆ Check that the image quality criteria described in the installation and maintenance manual are met before using the system.
- ◆ Maintenance work should be performed by trained technical personnel only.
- ◆ Purchasing a maintenance contract will keep the system in optimum condition.

This maintenance plan includes the required tests according to IEC 62353.

### Periodic maintenance

Periodic maintenance includes:

- Safety inspection
- Preventive maintenance
- Quality and function tests
- Replacement of safety-relevant parts subject to wear

Maintenance interval: 24 months for the entire system

The maintenance interval is estimated based on the workload of 40,000 patients within 24 months. If the examination workload is higher than the estimation, the maintenance interval may be shortened.

This work may be performed only by qualified and authorized service engineers. Qualified in this context means that the engineers have been trained accordingly or have acquired the necessary experience in practice. Authorized means that the engineers have been authorized by the operator of the system to perform maintenance work.

Prior to the first system startup, we recommend that you name a staff member who will be responsible for ensuring that routine checks and preventive inspection and maintenance work are performed. This employee has to archive all certificates.

In addition to our repair service, Siemens Healthineers also offers you a complete range of services for preventive inspection and maintenance of your system. These services can be called on as required or agreed in a flexibly drafted service contract.

If you have not received a quotation from our Customer Services organization, please contact the Siemens Healthineers representative responsible for your facility.

**Safety inspection** The following checks are intended to ensure the safety of the system. Where appropriate, preventive measures must be adopted or repairs performed. For the most part, the points to be checked are prescribed by laws and standards.

Maintenance interval: 24 months for the entire system

Listing of work steps to be performed:

| Object or function                  | Reason                                                                                     | What is checked                                                                                                                                                                                                                        |
|-------------------------------------|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>General checks</b>               |                                                                                            |                                                                                                                                                                                                                                        |
| Overall system                      | Safety of patients and personnel                                                           | Visual check of the system (cover panels) for damage and sharp edges                                                                                                                                                                   |
| Cables and how cables are laid      | Protection of patient and personnel against electric shocks                                | Visual check of cables and how cables are laid for safety deficiencies                                                                                                                                                                 |
| Accessories                         | Safety of patients                                                                         | Check for safety deficiencies                                                                                                                                                                                                          |
| System radiation protection devices | Protection of patient and personnel against radiation injuries                             | Check unit radiation safety devices as well as any additional safety devices that are configured, e.g. lower body radiation shield, upper body radiation shield, floor-mounted radiation shield, for installation and for damage.      |
| Operator documents                  | Protection against incorrect operation and thus protection of patients, personnel and unit | Check for the presence of operator documents                                                                                                                                                                                           |
| Warning notices                     | Protection against incorrect operation and thus protection of patients, personnel and unit | Check for presence of the required warning labels visible to the user for operation of the system.                                                                                                                                     |
| Insertable fuses                    | Protection of patients, personnel and unit                                                 | Check all insertable fuses that are accessible without tools to determine whether manufacturer data (nominal value and switch-off characteristics) are met.                                                                            |
| User interfaces                     | Protection against incorrect operation and thus protection of patients, personnel and unit | Check for legibility and operating symbols.                                                                                                                                                                                            |
| <b>Electrical checks</b>            |                                                                                            |                                                                                                                                                                                                                                        |
| Electrical safety                   | Protection of patient and personnel against electric shocks                                | <ul style="list-style-type: none"> <li>• Measurement of ground wire resistance</li> <li>• Device leakage current measurement</li> <li>• Patient leakage current measurement of the complete system per national regulations</li> </ul> |

| Object or function                     | Reason                                                                                     | What is checked                                                                                                                                                                                        |
|----------------------------------------|--------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Mechanical checks</b>               |                                                                                            |                                                                                                                                                                                                        |
| Floor mounting                         | Protection of patients, personnel and unit                                                 | Mounting (visual check)                                                                                                                                                                                |
| Pull and support cables                | Protection of patients, personnel and unit                                                 | Check for wear and fraying                                                                                                                                                                             |
| Chain drives                           | Ensure power transmission and safety functions                                             | Check for wear, tension, ease of movement and function of the guide rollers                                                                                                                            |
| System movements (manual)              | Safety of the system latching mechanisms                                                   | Check of brake holding and end stops                                                                                                                                                                   |
| Movable components                     | Protection of patients, personnel and unit                                                 | Cable lead-in, movement behavior and, if applicable, brakes                                                                                                                                            |
| Accessories                            | Safety of patients, personnel                                                              | Check for functionality and for good condition                                                                                                                                                         |
| User interfaces                        | Protection against incorrect operation and thus protection of patients, personnel and unit | Check for damage                                                                                                                                                                                       |
| Safety-relevant parts subject to wear  | Protection of patients, personnel and unit                                                 | Check for damage                                                                                                                                                                                       |
| <b>Function check</b>                  |                                                                                            |                                                                                                                                                                                                        |
| Emergency stop devices                 | Preventing the first malfunction caused by application errors and collapsing patients      | Switch-off of system functions after activating the emergency stop device                                                                                                                              |
| Control devices and warning indicators | Inform the operator about relevant system conditions and overload situations               | The functions of the following displays: <ul style="list-style-type: none"> <li>• Radiation</li> <li>• X-ray tube load</li> <li>• Blocking</li> <li>• Unit movements at the collision limit</li> </ul> |
| System movements (motorized)           | Protection of patients, personnel and unit                                                 | Safety shutdown of the movements                                                                                                                                                                       |
| Collision protection                   | Prevent damage to unit components                                                          | Automatic shutdown of system movements in the collision zone (for example ceiling, wall, floor)                                                                                                        |
| Function test                          | Prevent damage to unit components                                                          | Final function test of all components                                                                                                                                                                  |

\* DHHS and country-specific regulations must be observed

† Depending on the system configuration

**Preventive maintenance** The purpose of preventive maintenance is to reduce unforeseen failures to a minimum. Thus the prerequisites for the system to meet the assured characteristics in the long term are created.

The effects of different operating conditions (full or partial load operation, temperature, size of dust particles, humidity, gases, vapors) are checked and the condition of parts subject to wear is determined by recording and analyzing characteristic values. Preventive measures must be adopted or repairs must be undertaken as appropriate.

The specified maintenance intervals correspond to the minimum requirements. Compliance with stricter national legislation may be necessary.

Maintenance interval: 24 months for the entire system

Listing of work steps to be performed:

| Object or function | Reason                                                                                                                                                                    | What is checked                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
|--------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Overall system     | Preventive measures to avoid: <ul style="list-style-type: none"> <li>• Safety risks</li> <li>• Overheating</li> <li>• Wear and tear</li> <li>• Image artifacts</li> </ul> | <ul style="list-style-type: none"> <li>• Check of operating data</li> <li>• Check of cables, cable routing and cable connections for damage</li> <li>• Cleaning contrast medium, blood, disinfectants from areas not accessible to the customer</li> <li>• Inspection and cleaning of coolant and air circulation passageways</li> <li>• Inspection and cleaning of optical transmission paths</li> <li>• Removal of foreign objects, for example positioning aids and injection needles</li> <li>• Touch-up of paint to prevent corrosion and infection</li> <li>• Check of measurement values with tolerance ranges</li> <li>• Check of the movement forces</li> <li>• Check of the drive characteristics, acceleration and deceleration movements</li> <li>• Measures to ensure that all components move smoothly</li> <li>• Check and analysis of points of friction</li> <li>• Reading and analyzing the error log files</li> <li>• Repairing minor damage</li> </ul> |

#### Quality and function tests

Quality and function tests are used to check whether the system meets the assured characteristics. Image quality tests determine differences from the original status (for example resolution, contrast range, minimum contrast, image signal and, if applicable, check of Digital Subtraction Angiography).

If there are differences, preventive measures must be taken or repairs performed, whenever is appropriate.

Maintenance interval: 24 months for the entire system

Listing of work steps to be performed:

| Object or function                          | Reason                                                                                                                                                                                                                   | What is checked                                                                                                                                                                    |
|---------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Overall system                              | Optimum operation on the basis of the specifications listed in the data sheet                                                                                                                                            | Interaction of all system components according to the assured characteristics                                                                                                      |
| Power line connection                       | Unrestricted power input from the line power source for the purpose of maximum system load. (The power input can be limited by oxidation and corrosion, which can result in exposure fluctuations and system shutdowns). | Internal line impedance                                                                                                                                                            |
| Vacuum components:<br>• X-ray tube assembly | Assuring the system specifications (vacuum components are subject to aging)                                                                                                                                              | Image quality                                                                                                                                                                      |
| Beam geometry, centering, beam collimation  | Observance of the specifications and legal regulations* to minimize radiation exposure of patient and personnel                                                                                                          | Centering the central beam onto the center of the image receptor.<br><br>Coincidence of radiation field size and image receptor size or light field size and radiation field size. |
| Radiation dose                              | Compliance with specifications and legal requirements* to minimize radiation exposure of patient and personnel                                                                                                           | Check of dose rate/cutoff dose† during<br>• Acquisition in all operating modes                                                                                                     |
| Detail recognition                          | Compliance with specifications and legal requirements*<br>Ensuring image quality                                                                                                                                         | Resolution† during<br>• Acquisition in all operating modes                                                                                                                         |
| Image contrast                              | Compliance with specifications and legal requirements*<br>Ensuring image quality                                                                                                                                         | Minimum contrast and dynamic range† during<br>• Acquisition in all operating modes                                                                                                 |
| Image display                               | Compliance with specifications and legal requirements*<br>Ensuring image quality                                                                                                                                         | Brightness, focus and the geometry of the configured displays                                                                                                                      |
| Image artifacts                             | Compliance with specifications and legal requirements*<br>Ensuring image quality                                                                                                                                         | Image display for intolerable image artifacts in all existing application techniques                                                                                               |
| Dose measurement device†                    | Maintenance of specifications                                                                                                                                                                                            | Accuracy of the display                                                                                                                                                            |
| Image documentation systems                 | Compliance with specifications and legal requirements*<br>Ensuring image quality                                                                                                                                         | Gray scale reproduction, geometric imaging characteristics, resolution, optical density, artifacts                                                                                 |

\* DHHS and country-specific regulations must be observed

† Depending on the system configuration

#### Replacement of safety-relevant parts subject to wear

Safety-relevant parts which are subject to wear must be replaced in periodic intervals

## 3D ceiling stand

| Object or function | Reason                                                                  | What is replaced?                         | Interval |
|--------------------|-------------------------------------------------------------------------|-------------------------------------------|----------|
| 3D                 | Wear can cause springs to break and the support cable to split or fray. | Spring-loaded mechanism with steel cables | 8 years  |

## 2.1.9 Installation, repair or modifications

Modifications of or additions to the product must be made in accordance with the legal regulations and generally accepted engineering standards.

Siemens Healthineers does not accept responsibility for the safety features and for the reliability and performance of the equipment as the manufacturer, assembler, installer, or importer, in the following cases:

- Installation, equipment expansions, readjustments, modifications, or repairs are not carried out by persons authorized by us.
- Components affecting safe operation of the product are not replaced by original spare parts in the event of a malfunction.
- The electrical installation of the room does not meet the requirements of the VDE regulation or the corresponding national regulations.
- The product is not used in accordance with the Operator Manual.

### Technical documents

On request, we can give you technical documents for the product for a charge.

These documents do not imply authorization to carry out repairs.

We recommend that you obtain a report showing the nature and the extent of the work performed from the persons carrying out such work. The report should include all changes in rated parameters or operating ranges as well as the date, the name of the company and a signature.

We accept no responsibility for repairs performed without our express written approval.

### Installation

All information and installation requirements contained in the service documentation require strict adherence.

If you require additional information, please contact the qualified people specified by the manufacturer.

### Replacement

You can replace the following parts if necessary once the system has been switched off:

- Keyboard
- Mouse
- Barcode reader (optional)
- IR sensor
- Hand switch for radiation release (optional)

You can also replace the following parts as required before examinations:

- Grid
- Detector battery

### 2.1.10 Error messages and status indicators

The system displays warnings and error messages as well as information about current system status on the screen. It gives information on how to proceed in each situation.

For further details refer to (→ Page 58 *Ceiling-mounted tube support*), or to the Imaging System Operator Manual.

## 2.2 Personal safety measures

### 2.2.1 Information about unit movements

Because of the possible speed of movement of the system, it must be operated very carefully.

#### WARNING

System movements cause collision

##### **Risk of injury or system damage**

- ◆ Perform system movements only when it is certain that neither the operator, the patient, third parties nor pieces of equipment can be endangered by these movements.
- ◆ Remove any objects or accessories, for example injector or infusion stand, from the collision area.
- ◆ Make sure you are standing outside the dangerous area.
- ◆ Always watch the patient when performing system movements.

#### CAUTION

Touching the column during movement

##### **Risk of crushing**

- ◆ Do not place your hands to moving or rotating parts.
- ◆ Use the provided hand grips to move ceiling-suspended parts.
- ◆ Whenever this is not possible, pay special attention to all crush zones on the column especially those involving the hands or fingers.

Unit movements should be started only if neither a patient nor a third party may be endangered and if no object may hinder unit movements.

However, if an object (for example, during a down movement) was hindering a movement direction (for example seat) or the height of the detector unit or the X-ray tube was manually readjusted, the alignment control must be readjusted.

If this alignment is not performed, no acquisitions are possible!

## Deadman function

All unit movements are controlled with a deadman grip (DMG), that is, movements are performed only while the operating element is being actuated. In the event of a danger, stop the movement immediately by releasing the operating element.

## Danger of crushing

The patient and operating personnel must grip only the handles which are intended for the proper handling of the equipment or positioning of the patient. Where this is not possible, attention must be paid to possible risks of injury by crushing in the vicinity of moving parts.

- Pay special attention to the risks of crushing fingers/hands between moving parts and their guide openings.
- Before performing unit movements make sure that the patients do not grasp the frame of the tabletop.

## Abnormal movements

If any part of the system moves although movement was not released, there might be a fault.

Example:

The monitor suspension system moves down by itself.

- Shut down your system and call Siemens Healthineers Service.

## Collision protection

The MULTIX Impact C has a collision control device comprising all motor driven moving components like patient table and wall stand.

Attention: When organ programs are used in which the tube was programmed in the collision zone of the table, collisions of the tube and the table are possible.

Never leave the system unattended during system movements.

## Floating tabletop

Please pay attention to the behavior of the floating tabletop when the system is turned off.

### WARNING

Brakes on tabletop have reduced force when the system is turned off.

#### Risk of injury

- ◆ Check that nobody is leaning on the tabletop when the system is turned off.
- ◆ Avoid positioning patients such that a sudden movement of the tabletop could cause them to fall.
- ◆ Execute patient transfer with special care.

## Red emergency STOP buttons

### WARNING

Movement without intentional activation of control elements.

**Risk of collision, risk of injury to patient or operator, risk of damage to the system.**

- ◆ If movement does not stop, press the nearest emergency STOP button.

### Triggering STOP

If a malfunction of a unit movement causes an emergency situation, danger to the patient, to operating personnel, or to the unit, do the following:

- ◆ Press one of the red emergency STOP buttons immediately.

All system motor drives are shut down and movements are stopped immediately. Motorized movement can only be continued when STOP is canceled.

Radiation is not interrupted, and can be released again even without canceling STOP.

### Canceling STOP

Only when the cause of the danger has been unequivocally identified and remedied, can the emergency STOP button be disengaged.

- ◆ To do so, turn clockwise and pull the red emergency STOP button.


In case of a system failure, press the emergency STOP button and unlock it again. The system is reinitialized.


## Where are the emergency STOP buttons?

You find the emergency STOP buttons in the following locations:


On the table

On the generator ON/OFF console

On the wall stand (The emergency STOP buttons are located at both sides of the wall stand.)

### Emergency SHUTDOWN button (installed on-site)

Use this shutdown method only in an extreme emergency because it is an uncontrolled process!

Consequences:

- All unit movements are interrupted.
- Radiation is turned off.
- The current system program is interrupted.
- Current operating sequences are interrupted and cleared.
- All current acquisition data are cleared, unless they have been saved to non-volatile memory.

Data could be lost, for example unsaved images, exporting and filming jobs, and so on. The tube cooling system is also disconnected from the power supply, and the tube can overheat.

If a flat detector is used, it is also disconnected from the power supply. After the power has been switched on again, an additional waiting time is required to ensure optimum image quality.

An emergency power supply, if one is connected, is not activated when emergency SHUTDOWN is pressed.

The UPS of the imaging system is also not activated.

**Shutdown in an emergency** Only exclusively, if there is danger for patients, operators, third parties or the unit:

**CAUTION**

Shutdown using the emergency SHUTDOWN button

**Risk of data loss and/or damage to tube due to lack of cooling**

- ◆ Use the emergency SHUTDOWN button only in case of emergency or when the system cannot be switched off normally.

- ◆ Press the on-site emergency SHUTDOWN button.

The entire system is disconnected from the power supply.

**Switching on again**

Only if the cause of the danger has been unequivocally identified and remedied, the emergency SHUTDOWN button can be released and the system operated again.

- 1 In all other cases, for example system malfunction, contact Siemens Customer Service.
- 2 If the stop function does not respond normally, immediately activate an emergency shutdown to turn off the entire system.

If this happens, you may not continue to use the system. Contact the Siemens Customer Service.

**Emergency procedures****CAUTION**

Medical devices are not infallible and unexpected errors may occur

**System is not available in an emergency or may stop working during an examination**

- ◆ Have an emergency plan in case of non-availability of the system, for instance, to use a different system.
- ◆ Have an emergency plan for continuing or ending different types of examinations if the system stops working, for instance, during endoscopy, use the viewer of the endoscope to continue or end the examination.
- ◆ Consider whether an uninterruptible power supply is necessary at your site.

**Power failure**

If the power network at the site is very unstable, it may be advisable to install an uninterruptible power supply to prevent data loss. Contact Siemens for more information.

With UPS for imaging system (Uninterruptible Power Supply)

- 1 In case of a short power failure, please try to switch on the system after some seconds.
- 2 If it was possible to switch on the system, please wait until the system is ready.

- 3 Then continue working or rescue the patient.
- 4 If it was not possible to switch on the system or if the patient has collapsed, remove the patient immediately.

For more details on the UPS, refer to (→ Page 83 *UPS (for the imaging system)*)

## 2.2.2 Danger zones / danger points

### ► Danger points

The positions marked in the following illustrations show danger zones in which patients or operating personnel could suffer injury by crushing or colliding:

- Due to the mobility of the stands, there is a danger of injury in certain areas during unit movements.
  - Make sure that there are neither persons nor objects in the swivel or rotation range of the stands.
  - Make sure that there are neither persons nor objects in the lifting range of the patient table.
  - Make sure that the knees of sitting staff members are not in an area under a stand or the patient table.
  - Avoid standing or sitting immediately adjacent to the system.
  - If it is necessary to position the patient with legs or knees under the crossbeam at the head or foot end of the table, avoid system movements.
- Be aware of the danger zone between patient table and detector tray.

### WARNING

Unintended overtravel when lowering the table

**Crushing of patient or operator seated next to the table (for example, a wheelchair patient)**

- ◆ Avoid movements when somebody is within the dangerous movement zone.
- ◆ Table movements do not stop immediately when the control elements are released.
- ◆ Most axes perform a short overtravel movement.

## Patient table


Fig. right side: Patient sitting with knee under the patient table (not allowed)

Above figure on the right shows a particular patient positioning that is not allowed because of the potential hazards when lowering the patient table.

Please pay more attention to this warning sign placed on the foot end and head end of the table.


### WARNING



Lowering the table onto a seated patient

**Crushing of patient or operator seated with legs under the table (for example, a wheelchair patient)**

- ◆ Do not position a patient with legs under the table.
- ◆ Do not move the table while a patient is seated near the table. System will stop when collision is detected, but there is still a risk of injury for the patient!
- ◆ Move the patient away from the table first.


## X-ray tube support


**Wall stand****Wallstand with pulled out detector tray**

### Lateral detector holder

### Lateral cassette holder

### Grip locations

When handling the system correctly and when positioning the patient, operators and patients must use only the grip locations provided for this purpose.

The following grip locations are provided:

- Two hand grips (front and back)

#### Grip attachment

- ◆ Make sure that the hand grips are always attached.

If these grip locations cannot be used:

Pay special attention to the mentioned possibilities of crushing between moving parts and their guide openings.

Make sure during the examination that the patient does not hold on to the edges of the patient table under any circumstances.

## 2.2.3 Safety during patient examination

### Patient positioning

- All safety-related equipment must be installed and functioning, in particular the hand grips and the compression belt.
- The hands, arms, legs, head, and hair of the patient must not extend unsecured beyond the edge of the tabletop.

Observe the patient while moving the tabletop and during system movements.

#### WARNING

Patient cannot maintain position or becomes unconscious

##### Injuries to patient from falling

- ◆ Check that the patient is physically able to hold the examination position, especially when standing.
- ◆ Choose a different position or use positioning aids for immobilizing weak or very unstable patients.

#### WARNING

Hand grips loosen or are not used

##### Injuries to patient from falling

- ◆ Check that the patient hand grips on the table are in the right position.
- ◆ Ensure that hand grips are tightened prior to use.
- ◆ Make sure that the patient uses the hand grips during system movements or when turning over.
- ◆ Use other positioning aids for immobilizing weak or very unstable patients.

### Correct image orientation - patient orientation

It is the responsibility of the operator to ensure correct image orientation on the monitor or film.

**CAUTION**

Image is incorrectly flipped due to incorrect patient position information or right/left labeling.

**Risk of incorrect diagnosis or treatment because of left/right or up/down confusion**

- ◆ The examiner is responsible for using the image orientation functions and interpreting the images correctly.
- ◆ Make sure that the R/L labels have been placed correctly.
- ◆ Compare the patient position data with the anatomic view when diagnosing images to exclude any errors.
- ◆ If necessary, use lead letters or similar devices during acquisition.

## Visual contact to patient

Make sure that there is visual and acoustic contact with the patient. Thus, you stay informed about the condition of the patient at all times.

## Room lighting

According to DIN 68 68-157 (valid in Germany), the lighting in rooms in which diagnoses are made on image display devices (monitors) must fulfill the following requirements:

- The lighting must be adjustable and glare-free.
- The setting of the illuminance must be reproducible (for example, dimmer with scale).
- No mirroring or reflections of windows, lights, view boxes, and so on. must occur in the operating position of the monitors.

## Prerequisites for diagnosis and treatment planning

The MULTIX Impact C imaging system software has been designed and tested for use in diagnosis and treatment planning based on digital radiographic X-ray images and series.

To ensure that the imaging system produces monitor images suitable for these purposes, the monitor must meet certain criteria for image quality.

**Checks** Image quality can deteriorate over time because of aging and normal wear and damage of the monitor and other components.

- The image quality must be checked at regular intervals (once per month) after installation.
- Thus, you make sure that the system is still suitable for use in diagnosis and treatment planning.

For further information, refer to (→ Page 21 *Maintenance plan*).

The operator must ensure that qualified personnel are chosen and that the criteria for image quality specified in the installation and maintenance instructions are fulfilled.

**Test images** Test images for calibration of the monitor and to test the quality output of the laser camera, are stored in the system.

For further information, refer to the Imaging System Operator Manual.

## Displays

Please observe the following information:

- The operating indicator must light up.
- Please keep the ventilation slots of the monitors unobstructed at all times.

Do not touch the surface of the screen using sharp, pointed objects.

Do not open the display under any circumstances.

### CAUTION

The standard 24" display and 19" display (optional) are not intended for diagnosis.

#### Incorrect diagnosis possible

- ◆ For diagnosis use a diagnostic display or transfer the image to a workstation.

Monitors are suitable for medical on-line diagnosis only if special measures to assuring image quality are adopted (especially determining the brightness and contrast values).

In case of doubt, film images on a hardcopy camera if you want to make a diagnosis.

Siemens Healthineers undertakes no liability for diagnoses which are performed on non-Siemens Healthineers monitors or on monitors not calibrated by Siemens Healthineers.

## Use of hardcopy cameras

Only hardcopy cameras approved by Siemens Healthineers may be used. Approvals by Siemens Healthineers refer to DICOM image type XA. Siemens Healthineers does not accept any liability for diagnoses made on the basis of images from non-approved hardcopy cameras.

## Use of the full-field light localizer

### CAUTION

Photobiological effect of ultraviolet radiation

#### Eye injury

- ◆ Do not look into the light beam for longer than 15 seconds.
- ◆ Always keep enough distance to the collimator.

## 2.2.4 Radiation protection

The X-ray equipment MULTIX Impact C with radiation protection complies with IEC 60601-1-3 :2013, IEC 60601-2-54:2018.

Mode of operation: Continuous operation with intermittent loading

### Important information

The automatic format collimation system for acquisition helps reduce the radiation dose to the patient and examiner considerably.

For radiation dose indication during normal use of the equipment, the system is equipped with a measurement system based on an ionization chamber calibrated in dose area product.

#### Measures for reducing radiation exposure

Please observe the following:

Radiation may be released only by authorized operators or persons with proper authorization to apply ionized radiation.

Please pay attention to the fact that the air kerma is increased by the following settings:

- Higher kV
- Higher mAs
- Higher ms
- Smaller SID

#### Avoiding unwanted radiation

#### CAUTION

Tube is not directed toward the detector during acquisition.

#### Unwanted radiation exposure

- ◆ Use the free exposure mode with care.
- ◆ Check that the tube is active and use the light marker for positioning before releasing X-ray.

#### Radiation protection for the patient

- 1 If possible, ensure the best possible protection of the patient during acquisitions in the vicinity of the reproductive organs (use gonadal shields, lead-lined rubber covers).
- 2 Keep the radiation field as small as possible without reducing the active measuring field.
- 3 Remove all radiopaque parts from the radiographic field, if possible.
- 4 Set the X-ray tube voltage as high as possible (not forgetting the image quality; the optimum is 63 kV for iodine contrast).
- 5 Set the X-ray tube to skin distance as large as is reasonable for each examination.  
A larger focal spot distance leads to a lower air kerma.

**Radiation protection for the examiner**

- 1 If possible, release the acquisition series in the control room.
- 2 Stay in the examination area as short as possible.
- 3 If the exposure is to be released in the examination room, use an additional radiation shield (protection for upper and lower body) or radiation-proof window.  
These measures contribute greatly to your personal radiation protection.
- 4 If the protective shields are not used, wear radiation protection clothing with a 0.25 mm lead layer or similar clothing.
- 5 Keep the maximum distance from the source of radiation.
- 6 Check your personal dose by wearing a radiation monitoring badge or a pen dosimeter.

**CAUTION**

Leakage, extrafocal, stray, scatter and residual radiation

**Undesired radiation exposure**

- ◆ When releasing the x-ray, stand as far away as safely possible from the radiation beam path.
- ◆ Use radiation protection as necessary.

**CAUTION**

Failure to observe radiation protection regulations

**Unnecessary radiation exposure to personnel**

**This system produces x-rays.**

- ◆ Wear protective clothing such as lead aprons.
- ◆ Monitor radiation exposure using personal dosimeters.

**Switching off in an emergency**

If a problem, for example a defect, occurs during the examination and radiation can no longer be interrupted by letting go of the radiation release button:

- ◆ Press the nearest emergency STOP button.

## 2.3 Equipment safety measures

### 2.3.1 Protection against electric shock

To avoid the risk of electric shock, a protective conductor must be implemented when connecting this device to mains power.

## Power supply

For all products that are operated within an X-ray system, the power supply must be made through a contactor or other multipole circuit breaker installed on-site.

The room installation must comply with DIN VDE 0100-710 or the corresponding national regulations.

## Covers

If socket covers (especially socket covers of the operating modules) are damaged, they must be replaced.

In the event of defects, for example, if a covering cap has broken off, call the Siemens Healthineers Service.

## Protection class

The system belongs to protection class I with a type B applied part according to IEC 60601-1.

Protection against ingress of water:

- Foot switches: IPx8
- MAX wi-D: IP43 or higher
- Core XL: IPx3
- Rest of system: IPx0

## Potential equalization

Systems for which potential equalization is recommended must only be operated in medical facilities where supplemental potential equalization has been installed and tested according to the specifications in DIN 57107/VDE 0107/6.81 section 5 (Federal Republic of Germany) or the relevant local and federal regulations.

## Opening the units

Only authorized service personnel are permitted to open the units.

### 2.3.2 Electromagnetic compatibility (EMC)

#### CAUTION

Radio interference from the system affects nearby equipment

**Risk of malfunction of nearby equipment.**

- ◆ Reorient or relocate the equipment or system or shield its location.
- ◆ Test the operation of the equipment in its new location.

Refer to (→ Page 237 Notice concerning electromagnetic compatibility (EMC))

### 2.3.3 Safety-relevant parts subject to wear

This system contains no safety-relevant parts subject to wear.

### 2.3.4 Product service life

Our medical engineering products are designed for operation under average conditions in accordance with the Operator Manual for a useful life of 10 (ten) years. If the products are operated for longer than this time, additional checks and possibly repairs beyond the usual maintenance procedures may be necessary in order to ensure the functional integrity and safety of operation of these products. Please talk to us about these measures early enough!

### 2.3.5 Mechanical safety

#### Maximum weight

##### CAUTION

Too much weight on system components

##### **Risk of damage or patient injury**

- ◆ Observe the maximum allowable patient weight including the weight of the accessories.

The permitted patient weight on your tabletop is stated on the label on the tabletop and in the technical data.

It is important that the load is distributed uniformly over the tabletop. Otherwise there is a risk of material deformation and system malfunction. The patient weight includes any parts permanently or loosely connected with the patient, such as equipment, prostheses, implants, plaster casts.

#### Example of incorrect use

A patient with a maximum weight sitting at the end of a fully extended tabletop.

#### Maximum load

The component, for example tabletop or accessory, must not be loaded with more weight than indicated.

## Equipment damage

### Avoiding equipment damage

#### WARNING

System movements cause collision

##### **Risk of injury or system damage**

- ◆ Perform system movements only when it is certain that neither the operator, the patient, third parties nor pieces of equipment can be endangered by these movements.
- ◆ Remove any objects or accessories, for example injector or infusion stand, from the collision area.
- ◆ Make sure you are standing outside the dangerous area.
- ◆ Always watch the patient when performing system movements.

Check system and accessories daily for damaged or loose parts. If problems are found, have the parts repaired before resuming patient examinations.

- 1 Before starting system movements, especially lowering the table, make sure that the movement range is free of obstructions.
- 2 Remove chairs, steps, stands, waste containers, and similar objects from the movement area.

No collision monitoring!

- 3 Do not place any objects or consumable material on the cover of the table support, on the detector cover and on the longitudinal guides of the stand carriage.

Considerable forces which can damage these objects in the area of movement of the systems occur during movements of the tabletop.

- 4 Do not stand at any place on the covers of the table support.

The covers can be deformed.

Components located underneath them can be damaged and thus lead to operating disturbances.

#### WARNING

Movement without intentional activation of control elements.

##### **Risk of collision, risk of injury to patient or operator, risk of damage to the system.**

- ◆ If movement does not stop, press the nearest emergency STOP button.

### Damage to the tabletop

A damaged tabletop is a potential hazard to the patient!

When the patient tabletop collides with an obstacle such as a bed or instrument table, hairline cracks can result.

Call Siemens Healthineers Uptime Service and let the tabletop be checked immediately if it is possible that it has been damaged (for example accidental collision with the bed of the patient)!

## Damage to generator

 **CAUTION**

Objects on top of the generator cabinet block ventilation

**Generator becomes inoperable because it is too warm.**

- ◆ Do not store objects on top of the generator cabinet.

## Elevating table

When the tabletop meets any obstacle during downward movement of the table, the tabletop is lifted and movement stops. The system displays an error message.

In this case, only upward movement is possible.

**Warning signs**

Special danger points are marked on the unit with a warning sign:

## Danger of crushing

These warning signs indicate possible danger of crushing for the patient and examiner.

## Cardiopulmonary resuscitation

This warning sign shows the position of the patient table in cardiopulmonary resuscitation (CPR) with pressure compression up to 500N (50 kg).

## Risk of collision

Under certain circumstances, additionally installed components, for example lead protective screens, lamps, auxiliary equipment and so on, can cause collision with the system.

To protect these components and also the patient, such components are provided with the warning sign shown on the left.

We can provide you with these warning signs. The Siemens Healthineers Uptime Service will also attach them to additionally installed components on request.

## General danger

Caution, always observe the Operator Manual!

This warning sign indicates that special instructions can be found in the Operator Manual.

Please read the Operator Manual.

**Alert tones**

Alert tones will be played on the TUI when the system is in one of below situations:

- All components are ready for use after reboot.
- An exposure is released.
- X-ray tube is too hot (tube heat unit  $\geq 65\%$  HU).
- A confirmation messages is displayed on the TUI.
- The tube unit is approaching a detent.
- Top or center alignment is reached.
- The detector unit in BWS is moved to a predefined position via WRCC (optional).
- Collision is detected between tube unit and patient table.
- Motorized movement limit is reached.
- Skin dose threshold is reached ( $SID \leq 44$  cm).

### 2.3.6 Combination with other products/components

To ensure the necessary safety, only products or components expressly approved by Siemens, may be used in combination with the system.

The person connecting additional equipment to the medical device is considered to be the system configurer and is therefore responsible for ensuring that the current system configuration complies with the relevant standards (for example system standard IEC/EN 60601-ff and/or other applicable standards).

In case of doubt, please consult your local contact person.

### Interfaces

Accessories connected to the analog or digital interfaces must be certified according to the respective IEC standards.

Example:

IEC 60950-1 for Information Technology Equipment and IEC 60601-1 for Medical Equipment.

Furthermore all configurations shall comply with the valid version of the system standard IEC 60601-1.

Everybody who connects additional equipment to the signal input part configures a medical system. This person is thus responsible that the system complies with the requirements of the valid version of the system standard IEC 60601-1.

If in doubt, consult the Customer Services department or your local representative.

### 2.3.7 Disposal

#### CAUTION

Improper disposal of hazardous waste material can cause damage to people and the environment.

#### **Risk of injury and environmental damage**

- ◆ Dispose of waste material according to the national industry standards.
- ◆ Take account of local regulations governing the disposal of the product.

On disposing of the complete system or parts thereof, currently valid environmental legislation must be observed.

Examples of environmentally relevant components are:

- Accumulators and batteries
- Transformers
- Capacitors
- Flat panel monitors
- Phantoms

For details contact your local customer service representative or your Siemens Healthineers regional office.

NOTICE: System components which are hazardous to persons or the environment must be disposed of with care and in compliance with legally binding ordinances.

For all countries of the EU:

In the member states of the European Union (EU) Siemens Healthineers will take back and will dispose of the packing material of your system.

Products bearing this symbol are subject to EC directive 2002/96/EC on waste electrical and electronic equipment (WEEE), amended by 2003/108/EC. This directive determines the framework for the return and recycling of used equipment as applicable throughout the EU.

For details about return and disposal of the product or its components or accessories please contact your local Siemens Healthineers customer service or your Siemens Healthineers regional office.

### 2.3.8 Cleaning and disinfection

- (→ Page 48 *Important general information*)
- (→ Page 49 *Cleaning and disinfection process*)
- (→ Page 50 *Special cleaning information*)
- (→ Page 52 *Disinfection*)
- (→ Page 52 *Sterilization*)

## Important general information

### Safety hints

Switch off the system and disconnect it from the main power prior to cleaning or disinfecting.

#### CAUTION

Use of harsh cleaning agents, liquids, or sprays.

##### **Risk of electrical hazard or damage to the system**

- ◆ Use only substances for cleaning and disinfection, which are recommended.
- ◆ Do not let cleaning liquids seep into the openings of the system (for example, air openings, gaps between covers).
- ◆ Observe the cleaning and disinfection instructions.

#### CAUTION

Inadequate cleaning and disinfection

##### **Risk of infection**

- ◆ Clean and disinfect all contaminated surfaces and all parts which may come (or have come) into contact with patient after each examination.
- ◆ Use the recommended cleaning agents and disinfectants.
- ◆ Protect the portable detector with a single-use plastic bag.

### Ventilation slots

- ◆ Keep the ventilation slots of all components unobstructed.

#### CAUTION

Objects on top of the generator cabinet block ventilation

**Generator becomes inoperable because it is too warm.**

- ◆ Do not store objects on top of the generator cabinet.

### Generator cabinet

#### CAUTION

Cleaning liquids drip into the generator cabinet.

**Electrical parts are damaged.**

- ◆ Make sure that no liquids are used at the top of the generator cabinet.
- ◆ Clean the generator cabinet covers very carefully, making sure that no liquid seeps into the cabinet.

### Detector protection

Avoid contact of detector with iodine-containing or other strongly dyeing substances.  
We recommend to use a protective cover during examination.

### Cleaning in general

#### Preparations

- ◆ Clean all contaminated parts and all parts that may come or have come into contact with the patient, directly or indirectly.

When cleaning and disinfecting, use suitable gloves for your own protection.

#### Equipment

- 1 Only use water or a lukewarm diluted household cleaning agent solution.

The use of other than the compatible cleaning agents can result in damage to the equipment.

- 2 Wipe system parts with a damp cloth (squeeze out wet cloth before using it) until all contaminations are removed.
- 3 Remove any watery residues immediately.

### Disinfection in general

Disinfection by means of vaporization is not allowed!  
Disinfection by means of UVC light is not allowed!

### Cleaning and disinfection process

Preparation for decontamination: Protective equipment for cleaning personnel (gloves, water repellent protective apron, face protection mask or protective glasses and mask)

Always adhere to the disinfectant manufacturer's instructions for use regarding handling and exposure time of the disinfectant solution.

Always work from top to bottom (top down) and from clean to dirty areas.

Use wipes only as long as they leave a continuous liquid film on the surface.

If blood is to be wiped off, discard the wipe and use a new one to continue the cleaning and disinfection procedure.

A non-protein-fixing disinfectant with proven cleaning effect based on peroxide compounds, certified by corresponding national authorities (e.g. EPA or VAH), is recommended (Dismozon plus, Hartmann #981257, active ingredients: Magnesium monoperoxyphthalate hexahydrate 95.8g/100g was validated), tap water/flowing water (20+/- 2°C, minimum drinking water quality), lint free wipe (Braun Wipes Eco #19726 was validated).

- 1 Prepare cleaning and disinfection solution according to the manufacturer's instructions.
- 2 Wipe the surfaces of the product thoroughly with sufficiently saturated wipes for cleaning.
- 3 Make sure that all surfaces, existing grooves and protrusions of the product can be reached during the cleaning process and are completely moistened.
- 4 Soak a wipe with tap water and wipe all surfaces of the product thoroughly.
- 5 Check for cleanliness; if soiling is still visible, repeat the above steps.
- 6 Thoroughly wipe the surfaces of the product one more time with a sufficiently soaked wipe for disinfection.
- 7 Make sure that all surfaces, existing grooves and protrusions of the product can be reached during the disinfection process and are completely moistened.
- 8 Follow the manufacturer's instructions regarding the exposure time of the disinfectant solution.
- 9 Let the system dry.

### Special cleaning information

Never use abrasive cleaners or cleaning agents with solvents (e.g. cleaning solutions, alcohol or spot removers), since they may damage housing surfaces.

#### DAP chamber (Option)

High voltage is present in the DAP ionization chamber. It is therefore very sensitive.

- ◆ Use a soft cloth to avoid scratching the chamber.

#### Mobile detector

#### WARNING

The battery of the portable detector is not resistant to fluids.

#### Electrical shock, damage to the detector

- ◆ Remove the battery from the portable detector before cleaning the detector.
- ◆ Check that the detector is waterproof before immersing it (symbol on detector ).
- ◆ The battery can be wiped off but it is not waterproof.

The charge contact plate on the back of the detector must be completely dry before the detector is returned to the docking station, in the Bucky tray of the table or the Bucky wall stand. Otherwise, the detector will not charge sufficiently and the charge contact surfaces may corrode.

It is recommended to clean the back of the detector first so that it has more time to dry while the front is being cleaned.

The following agents can be used to clean the Core XL:

- Cleanisept
- Cavi Wipes XL
- Sani Cloth AF3 (Wipes)
- Neutral Disinfectant Cleaner
- Quaternary Disinfectant Cleaner
- Super Sani-Cloth
- Sani-Cloth Plus
- Apesin AP100
- Dismozon plus
- Clorox Healthcare Hydrogen Peroxide (Wipes)
- Oxivir Tb Wipes
- Optim 1 Wipes
- Cavi Wipes Bleach
- Peridox RTU
- Sani-Cloth Bleach
- Clorox Healthcare Fuzion Cleaner Disinfectant
- Clorox Healthcare Bleach Germicidal Wipes

**Touch User Interface (TUI)** You must clean the TUI regularly because it becomes dirty with fingerprints:

- 1 Do not spray the TUI directly.
- 2 Wipe the TUI only with a clean, damp cloth.

**Monitors, touchscreen and keyboard** When cleaning the monitor screens or displays use only a moist cloth without any cleaning agents.

- 1 Wipe the screen.
- 2 Dry it with a soft cotton cloth.
- 3 Remove immediately any contamination, for example, contrast medium stains.

**Liquid-crystal displays (LCD)** Liquid-crystal displays (LCD) are very sensitive to mechanical damage.

- 1 Avoid scratches and shocks.

- 2 Remove drops of water immediately, longer contact with water discolors the surface.

**Camera window** Avoid using organic acids to clean the camera window.

**Parts mounted above the patient** ♦ Remove dust on parts mounted above the patient regularly.

**Accessory parts** ♦ Please note that for some accessory parts, there are special instructions on cleaning in the corresponding chapters.

For more information, see the System Operator Manual - Chapter Accessories and Auxiliary Devices.

Unless special instructions are given there, the equipment cleaning instructions are valid.

## Disinfection

**Surfactant disinfectants** For surface disinfecting of the equipment, we recommend using liquid solutions of a commercially available surfactant disinfectant based on the following compatible ingredient classes:

- Peroxide compounds

As a general rule, disinfectant sprays must not be used since the spray mist may penetrate into the equipment which means the safety of the equipment can no longer be ensured (damage to electronic components, formation of flammable air solvent vapor mixtures).

Only use liquid disinfectants and wipe with a moist wipe.

Some substances contained in disinfectants (for example, aldehyde) are known to be hazardous to health. Their concentration in the air must not exceed the statutorily defined limits.

Follow the manufacturers' instructions for use of these products.

## Sterilization

The system does not require sterilization.

## 3 System Overview

### 3.1 Product description

The MULTIX Impact C system is equipped with following components:

- Ceiling-mounted tube support with ACSS (Automatic Cassette Size Sensing) and Tracking function
- Bucky wall stand
- Image Station

Further components can be:

- Patient table
- A total of 3 different detectors are available:
  - MAX wi-D (wireless mobile detector with handle)
  - Core XL (wireless mobile detector with handle)
  - Core static (integrated detector)

The MAX wi-D as well as the Core XL can be used in the patient table, Bucky wall stand, and for free exposures.

The Core static can only be equipped with the Bucky wall stand.

Scatter radiation grid:

- The patient table and the wall stand are equipped with a stationary scatter radiation grid.
- For free exposures the MAX wi-D and Core XL can be equipped with a clip-on grid.

#### 3.1.1 Application capabilities

A MULTIX Impact C system may be used for projection of radiographic examinations of both adults and children (removable scattered radiation grid).

- On the patient table, X-rays can be taken on the recumbent or sitting patient from head to foot. Exposures in the region of the skull, of the spine, of thorax, lungs and abdomen as well as of the extremities are possible.
- On the wall stand, X-rays can be taken on the standing or sitting patient (chest, and so on).
- In addition, exposures onto MAX wi-D and Core XL in the form of free exposures are possible.


### 3.1.2 System configurations

The following system configurations are available:

- MULTIX Impact C with ceiling-mounted tube support (ACSS & Tracking):
  - Bucky wall stand with Core static
  - Bucky wall stand with Core static, and optional Core XL
  - Bucky wall stand with Core static, and optional MAX wi-D
  - Bucky wall stand with Core XL
  - Bucky wall stand with MAX wi-D
  - Bucky wall stand with Core static and patient table with Core XL
  - Bucky wall stand with Core static and patient table with MAX wi-D
  - Bucky wall stand and patient table to be used with one Core XL
  - Bucky wall stand and patient table to be used with one MAX wi-D


## 3.2 Basic system

### 3.2.1 System overview in the examination room


- (1) MULTIX Impact C ceiling support with control panel, tube assembly, TUI (Touch User Interface) and collimator
- (2) Patient table (optional)
- (3) Bucky wall stand

### 3.2.2 System overview in the control room

- (1) Image Station  
(Touch screen as an option is available.)
- (2) Keyboard and mouse
- (3) Generator ON/OFF console with radiation release button
- (4) Hand switch for radiation release (optional)
- (5) Wireless Remote Control Console (WRCC) (optional)
- (6) Barcode scanner (optional)

### 3.2.3 Image Station

The Image Station is a digital imaging system with main purpose to manage patients (studies), acquire, display, postprocess, film, and archive X-ray images.

The main component of the Image Station is an all-in-one computer with the Windows 10® operating system.

Do not touch the surface of the screen using sharp, pointed objects.  
Do not open the display under any circumstances.


Operate the touch screen (optional) with one dry and clean finger.

For further information refer to the Imaging System Operator Manual.

## 3.2.4 Generator ON/OFF console

- (1) LED indicator  
Lights up when power is ON  
Lights off when power is OFF
- (2) Generator ON button
- (3) Generator OFF button
- (4) Radiation release button
- (5) Radiation indicator  
Lights up when radiation is ON  
Lights off when radiation is OFF
- (6) Emergency STOP button
- (7) Hand switch (optional)
- (8) Detector calibration switch (at the back side of the console)

The exposure is released either with the release button (4) integrated in the control console or with the hand switch (7) with spiral cable.

Before releasing a radiation, check the selected acquisition data on the imaging system.

The selected acquisition data should be checked during exposure parameter settings.

### CAUTION

Wrong radiation parameters

#### Risk of undesired X-ray exposure due to wrong parameters

- ◆ Only suitably qualified personnel should use the system.
- ◆ Check the radiation parameters before releasing X-ray.
- ◆ Pay special attention to the exposure index (EXI).

## Release exposure with hand switch

### Preparations:

- ◆ Press the release button down to the 1st pressure point and hold it until the radiation release (approx. 1 s to 2 s).

The rotating anode is accelerated to nominal speed.

The required filament current is set.

After warming up the X-ray tube assembly and starting up the rotating anode, it is possible to release the exposure specifically.

### Radiation release:

- 1 Press the release button down fully and keep it pressed until the exposure is ended.

The radiation indicator on the control console lights up and an acoustic signal sounds during exposure release.

Letting the release button go ends the exposure immediately and the exposure can be underexposed.

- 2 Observe the following to increase the life of the X-ray tube assembly:

If possible always select 80% tube utilization at the generator.

In the programs with small focus and a tube voltage less than 70 kV, absolutely program 80% tube utilization.

Under 70 kV work if possible with the large focus.

Keep the preparation time short and press the exposure release button

Especially if you are working on workstations with a high patient throughput (emergency workstations), it is recommended to pay attention to the above notes to prevent premature wear and tear of the tube assemblies.

## 3.2.5 Unit overview

- Patient table  
(→ Page 72 *Patient table*)
- Ceiling-mounted tube support  
(→ Page 58 *Ceiling-mounted tube support*)
- Wall stand  
(→ Page 75 *Wall stand*)

### 3.3 Ceiling-mounted tube support

#### CAUTION

Heavy items (such as lead aprons) hung on ceiling-mounted devices.

**Parts could be damaged or fall down**

- ◆ Do not place any additional load on the ceiling-mounted devices.


#### WARNING

Damaged cables inside the lifting column of the overhead tube or detector.

**Overhead tube or detector can fall down. Risk of patient or user injury.**

- ◆ If the overhead tube or detector cannot be moved vertically, do not force the movement.
- ◆ Call Siemens Healthineers Service.

The 3D overhead support consists of the following components:


- (1) Lifting column with LED indicator (optional)  
LED green: System is ready for exposure.  
LED yellow: Exposures are being performed.  
LED blue: System is not ready for exposure.
- (2) X-ray tube unit
- (3) Touch User Interface (TUI)
- (4) Handles with buttons  
Use the handles if you want to bring the tube unit manually into another position (vertical, longitudinal, transverse or rotation).
- (5) Multileaf collimator (ACSS)
- (6) Collimator control panel

You can control the movements of the 3D ceiling mounted tube support with the aid of the handles and the touchscreen.

### 3.3.1 Touch User Interface (TUI)


The TUI is a color touch screen user interface for control of multiple functions. The touch screen corrects its display orientation from portrait to landscape, when the X-ray tube RHA angle is  $> 50^\circ$ .

The display on the TUI informs you about the position of the support (SID, tube angle), the generator settings, patient name, organ programs, and so on.

Example of landscape orientation (left) and portrait orientation (right)

#### Layout

#### Example

- (1) Patient information field
- (2) Control field
- (3) Movement button indicators
- (4) System status and information field

Press the TUI with one finger.

The icons/buttons shown in this chapter only outline their performance and should not be taken as a detailed "pixel-by-pixel" design.

### Buttons for movements and indications on the TUI

Next to the handles on the tube unit, there are four buttons on each side of the TUI.

The functions of the buttons are indicated by symbols and colors on the TUI.

| Vertical alignment of the tube unit |                                            | Horizontal alignment of the tube unit |  |
|-------------------------------------|--------------------------------------------|---------------------------------------|--|
|                                     |                                            |                                       |  |
| Button symbol                       | Function of the button                     |                                       |  |
|                                     | Tube rotation about the vertical axis      |                                       |  |
|                                     | Floating movement up and down              |                                       |  |
|                                     | Motorized movement up                      |                                       |  |
|                                     | Motorized movement down                    |                                       |  |
|                                     | Longitudinal movement                      |                                       |  |
|                                     | Longitudinal movement is in lock position. |                                       |  |
|                                     | Transverse movement                        |                                       |  |
|                                     | Transverse movement is in lock position.   |                                       |  |

| Button symbol                                                                     | Function of the button                       |
|-----------------------------------------------------------------------------------|----------------------------------------------|
|  | Tube rotation about the horizontal axis      |
|  | Centering tube to the detector in bucky tray |

### Patient information field

**Emergency**

202007201155151

(1)

(2)

(1) Patient information

(2) Help button

Before starting acquisition, make sure that the name of the patient on the TUI agrees with the name of the patient to be examined.

**Displaying help** You can display context-sensitive help information on the TUI.

1 Tap the ? icon.

The ? icon is highlighted.

The user interface is in help mode and content sensitive help information is available.

2 Pressing an element and the help information is displayed.

3 To exit the help mode, tap the ? again.

### Control field

You can select the organ programs and check/change the acquisition parameters.

|              |             |                 |                 |   |
|--------------|-------------|-----------------|-----------------|---|
| <b>Basic</b> | <b>More</b> | <b>Position</b> | <b>Chest PA</b> | ▼ |
| (1)          | (2)         | (3)             | (4)             |   |

(1) **Basic** tab

(2) **More** tab

(3) **Position** tab

(4) **Organ** program selection

There are several basic types of buttons on the TUI.

| Abbreviation | Type              | Explanation                                       |
|--------------|-------------------|---------------------------------------------------|
| P            | Push button       | Starts an action<br>The action terminates itself. |
| T            | Toggle button     | Toggles between ON and OFF                        |
| S            | Switch button     | Switches among different states                   |
| Pop          | Popup menu button | Opens a popup menu                                |

### Patient size and workplace selection

You can select the patient size and workplace by using the pop up menu.

**Patient size** The selected patient size is displayed with an icon.



| Icon                                                                                | Short description   |
|-------------------------------------------------------------------------------------|---------------------|
|   | Small/thin patient  |
|  | Medium patient      |
|  | Large/thick patient |

**Workplace selection** The selected workplace is displayed with an icon.


| Icon                                                                                | Short description |
|-------------------------------------------------------------------------------------|-------------------|
|  | Patient table     |
|  | Wall stand        |
|  | Free exposure     |
|  | Cassette exposure |

## Direct accessible functions

## SID tracking



| Icon                                                                                       | Short description                                                                                                                      |
|--------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------|
|           | Disabled: Preconditions are not fully fulfilled.<br>(→ Page 98 <i>SID tracking</i> )                                                   |
|           | Not pressed: inactive                                                                                                                  |
|           | Pressed: active and tracking on target                                                                                                 |
| <br>(red) | Active but tracking target lost<br>- or -<br>During tracking-movement<br>Tracking function can be deactivated in "lost" state as well. |






## Tube/detector tracking

| Icon                                                                                         | Short description                                                                                |
|----------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------|
|           | Not pressed: inactive                                                                            |
|           | Pressed: active and tracking on target                                                           |
| <br>(red) | Active but tracking target lost<br>Tracking function can be deactivated in "lost" state as well. |





## Top alignment

| Icon                                                                                | Short description                                                                           |
|-------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|
|  | Disabled: Preconditions are not fully fulfilled.                                            |
|  | Not pressed: center-aligned<br>Press the centering button to activate the center alignment. |

|                                                                                    |                                                                                     |                                                                                                                                                                                                                                                                                                                                               |
|------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                                    | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |    | <p>Pressed: top-aligned</p> <p>Press the centering button to activate the top alignment.</p> <p>Preconditions:</p> <ul style="list-style-type: none"> <li>• Wall mode is selected.</li> <li>• Tube RHA angle is <math>-90^\circ \pm 45^\circ</math> or <math>90^\circ \pm 45^\circ</math> depending on configured position of BWS.</li> </ul> |
| Restore OGP parameter                                                              | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |    | Restore the organ program parameters to default settings.                                                                                                                                                                                                                                                                                     |
| Control element disable function                                                   | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |   | Not pressed: inactive                                                                                                                                                                                                                                                                                                                         |
|                                                                                    |  | <p>Pressed: active</p> <p>All buttons for motorized movement of the tube stand and other components such as BWS and patient table are disabled!</p>                                                                                                                                                                                           |
| Software stop (only available for systems with Ortho or SmartPositioning function) | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |  | Disabled                                                                                                                                                                                                                                                                                                                                      |
|                                                                                    |  | <p>Pressed: active</p> <p>Software stops at <math>-90^\circ</math>; <math>0^\circ</math>; <math>+90^\circ</math> are turned on for tube unit rotation around the horizontal axis.</p>                                                                                                                                                         |
| Tube position                                                                      | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |  | Tube rotation horizontally (RHA)                                                                                                                                                                                                                                                                                                              |
| Generator functions                                                                | Icon                                                                                | Short description                                                                                                                                                                                                                                                                                                                             |
|                                                                                    |  | kV value selection                                                                                                                                                                                                                                                                                                                            |

| Icon                                                                               | Short description    |
|------------------------------------------------------------------------------------|----------------------|
|   | ms value selection   |
|   | mAs value selection  |
|   | Dose value selection |
|   | Density correction   |
|   | AEC fields selection |
|  | Focus selection      |

## Collimator functions

| Icon                                                                                | Short description                                                                |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------------------|
|  | Cu filter selection                                                              |
|  | Restore collimator button<br>(→ Page 138 <i>Restoring the last collimation</i> ) |
|  | Collimator size selection                                                        |
|  | SID adjustment for free exposures                                                |

## System status and information field

## System status

| Icon                                                                                | Short description                                                    |
|-------------------------------------------------------------------------------------|----------------------------------------------------------------------|
|  | Status: Ready                                                        |
|  | Status: not Ready<br>Error messages display with the warning symbol. |

| Grid status                                                                                             | Icon                                                                                                                                                                                                                                                                                                         | Short description |
|---------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-------------------|
| <br>Current: XXX       | A grid is inserted in the patient table or wall stand, and the grid type is indicated.<br>Grid type: <ul style="list-style-type: none"><li>• Universal: universal grid, <math>f_0= 140</math> cm</li><li>• 115 cm: grid, <math>f_0= 115</math> cm</li><li>• 180 cm: grid, <math>f_0= 180</math> cm</li></ul> |                   |
| <br>Current: none      | No grid is inserted.                                                                                                                                                                                                                                                                                         |                   |
| <br>Current: unknown   | For free exposures, a clip-on grid may be used.<br>The status of the clip-on grid cannot be detected by the system.                                                                                                                                                                                          |                   |
| <br>Current: removed   | Grid removed during grid test<br>Table or wall mode is active, grid runs out of its position during grid test.                                                                                                                                                                                               |                   |
| <br>Desired:XXX      | The grid defined in the organ program is indicated.                                                                                                                                                                                                                                                          |                   |
| <br>Desired: none    | No grid is defined in the organ program.                                                                                                                                                                                                                                                                     |                   |
| <br>Desired: unknown | An organ program for free exposure has been selected.<br>The grid is not defined in the free exposure organ programs.                                                                                                                                                                                        |                   |

The grid information is highlighted in yellow when the actual grid status does not match the setting in the organ program.

| Tube position | Icon                                                                                | Short description       |
|---------------|-------------------------------------------------------------------------------------|-------------------------|
|               |  | SID in cm or inch       |
|               | <b>RHA: 0°</b>                                                                      | Tube rotation angle RHA |

## Centering display

| Icon                                                                              | Short description                            |
|-----------------------------------------------------------------------------------|----------------------------------------------|
|  | Beam in center- or top aligned mode          |
|  | Red: beam not in center- or top aligned mode |

## Collimation (for automatic collimator only)


| Icon                                                                                   | Short description                                                  |
|----------------------------------------------------------------------------------------|--------------------------------------------------------------------|
|  30cm | Collimator height in cm or inch                                    |
|  30cm | Collimator width in cm or inch                                     |
|       | Collimator rotation<br>It displays when the collimator is rotated. |

## Collimator mode (for automatic collimator only)


| Icon          | Short description |
|---------------|-------------------|
| <b>ACSS</b>   | ACSS mode         |
| <b>Manual</b> | Manual mode       |

## 3.3.2 Multileaf collimator

## Control elements at the front

Collimator AL04 II-D

#### Collimator RFU

- (1) Adjusting knob for format width collimation  
Turning the knob to the left or right
- (2) Adjusting knob for format height collimation  
Turning the knob to the left or right
- (3) Button for full-field light localizer and laser line light localizer  
Switches off automatically after a preset time (configurable)
- (4) Tape measure for SID setting (cm and inch)
- (5) Accessory rails with locking spring
- (6) Detent lever for  $\pm 45^\circ$  rotation of the collimator about the vertical axis. The collimator only stops in the  $0^\circ$  position.


#### Additional Cu filter

There are four possibilities of selecting internal additional Cu filters for the automatic collimator:

- 0.0 mm (no) Cu
- 0.1 mm Cu
- 0.2 mm Cu
- 0.3 mm Cu

The selection of additional filters is specified in the organ program. However, you can change this setting at the imaging station or on the TUI.

#### Selecting additional filters

- ◆ Touch this button on the **Basic** tab on the TUI.  
A popup menu appears where you can select the required filter

The additional filter which is selected is displayed at the imaging station.

When changing the additional filters with respect to the organ program selected, this is marked with a star \* at the imaging station, denoting the change.

– or –

Select the filter at the imaging system.

Refer to the Imaging System Operator Manual for further information.

**CAUTION**

Incorrect filter selection

**Risk of increased radiation dose for the patient**

- ◆ Select the filter carefully.

## Accessory rails

The accessory rails may be loaded with maximum 70N (7kp).

**WARNING**

Falling objects because too much weight is hung on the collimator

**The rails could be damaged by the use of excessive force. A proper seating of the accessories or auxiliary devices can no longer be guaranteed and can lead to patient injury. The position of the tube assembly can change.**

- ◆ Do not hang more than 7 kg on the collimator accessory rails.
- ◆ Do not exert any force on the multileaf collimator or its rails when inserting the accessory or auxiliary device.
- ◆ Check that the accessory is firmly anchored in the collimator rails.


## Measuring tape for SID/SOD measurement

The measuring tape (6) has two scales:

Inches (left scale) and cm (right scale).


### Control elements at the underside

- (1) Laser line light localizer with slider to cover the outlet opening
- (2) Crosshairs in the light localizer window
- (3) Locking spring for accessories
- (4) Camera (optional)

#### Laser line light localizer

The laser line light localizer projects the axis mark for longitudinal centering, which is aligned with the centering mark on the receptor.

The laser line light localizer is switched on and off jointly with the full field light localizer.

Do not use optical lenses, mirrors and similar instruments if working with laser light. The optical instruments within the laser beam may amplify the laser intensity to dangerous values for eyes and skin.

Switch off the laser light when optical instruments are used.

To protect the eyes of the patient or any other person, the laser radiation outlet (1) of the laser line light localizer can be closed with the sliding cover (1a).

**CAUTION**

Looking into the laser light

**Eye injury**

- ◆ Take care that neither you nor anybody else looks directly into the laser beam.
- ◆ Close the sliding cover over the laser light to protect the eyes of the patient and other people.
- ◆ Note that the collimator can heat up when the light localizer is used for a longer period of time.

**CAUTION**

If the LED of the light localizer burns for a long time, parts of the light localizer can heat up.

**Danger of burns**

- ◆ Avoid contact with the light localizer to avoid burns.

### Switching on/off the laser line light localizer

- 1 To switch on the laser line light localizer, press the button at the front of the collimator.
- 2 To switch off the laser line light localizer, press the button again.

The laser line light localizer can also be switched off automatically by an internal time switch. (The internal time switch can be set by Siemens Healthineers service. 5 s to 118 s)

### Crosshairs

The crosshairs project the longitudinal and transverse axis of the radiation field onto the cassette or directly onto the patient.

- 1 The full field light localizer for projecting the crosshairs is switched on with the button at the front of the collimator.
- 2 To switch it off press the button again.

The full field light localizer can also be switched off automatically by an internal time switch. (The internal time switch can be set by Siemens Healthineers service.)

The laser line light localizer and the full field light localizer cannot be switched independently of one another.

### Locking spring

The locking spring is located on the left guide rail at the underside of the collimator.

The locking-spring locks the compensating filters, templates and so on, inserted in the accessory rails of the collimator, thus securing them against falling out.

### Removing the accessories

- ◆ Press the locking spring to the left until the compensating filter, template and so on, can be removed from the collimator.

## 3.4 Patient table

### 3.4.1 Application


The patient table is designed for use in stationary radiographic systems. The patient table supports the sitting or recumbent patient during the radiographic exposure and contains the detector tray. The tabletop is movable in the direction of the X-, Y- and Z-axis.

The patient table may be used for projection of radiographic examinations of both adults and children.

- X-rays can be taken on the lying or sitting patient
- Exposures in the region of the chest, of the spine, of thorax, lungs and abdomen as well as of the extremities are possible.

### 3.4.2 Overview

The patient table can be used with a MAX wi-D or a Core XL.


- (1) Tabletop
  - (2) Accessory rails
  - (3) Detector tray, movable to the left or right manually
  - (4) Emergency STOP push button on both front and rear sides
  - (5) Brake knob to move the detector tray
  - (6) Foot kick switch for tabletop movements and table lift
- As an option, a second set of foot kick switches on the rear side is available

### 3.4.3 Brake knob

The brake knob is mounted to the right side of the detector tray.

With the brake knob, the detector tray can be moved manually in longitudinal direction. The detector tray can be stopped and fixed at any position within its moving range.

### Moving the detector tray


- 1 Pull the brake knob (5) to release the brake of the detector tray.
  - 2 Let go of the brake knob in the desired position.
- The brakes now hold the detector tray in position.

The brake knob is not functional during the coupling of the tube unit and the detector unit in the fixed table.

#### 3.4.4 Detector tray with wireless detector

Example, detector tray with MAX wi-D

- (1) Wireless detector  
(→ Page 77 Wireless detectors (MAX wi-D or Core XL))
- (2) Handle of the detector tray
- (3) Stops

### 3.4.5 Emergency STOP push button


It interrupts all motor driven system movements.

For more details, please refer to (→ Page 14 System Safety).

### 3.4.6 Foot kick switches in the table base

#### Standard functions

- (1) Release the tabletop brakes (in longitudinal and transverse direction),
- (2) Table downwards movement, "Push button mode"
- (3) Table upwards movement, "Push button mode"

#### Alternative operation modes switch

The following modes of operation can be configured for both switches (1):

- Push button mode (dead man function)

The foot kick switch operates like a push button.

The corresponding function is enabled only as long as the foot kick switch is actuated.

- Sensor mode (toggle function)

The foot kick switch operates like a switch.

If the foot kick switch is actuated once, then the function is switched on (brakes released).

The function is switched off again at the second actuation.

#### **CAUTION**

Floating tabletop in sensor mode (toggle function) is not locked.

#### **Risk of falling or sliding**

- ◆ In case the floating tabletop is working in sensor mode, always lock the tabletop after patient positioning is finished.

## 3.5 Wall stand

### 3.5.1 Description

- The wall stand contains a detector unit with grid.
- The detector unit is fully counterbalanced and features continuous vertical adjustment.
- The detector unit can be continuously tilted from the vertical into the horizontal position (+90°) and up to -20° with a detent at 0° and 90°.
- Automatic synchronization (detector height adjustment)  
X-ray tube follows synchronized and keeps centering on middle of detector or its upper edge
- The wall stand has a centering function  
Fully automatic centering of X-ray tube on middle of detector or its upper edge
- The wall stand can be mounted on the floor.


### 3.5.2 Application

The wall stand can be used as examination unit for universal medical applications performed on standing, sitting or recumbent patients and involving exposures of the thorax, abdominal, pelvic, skull and vertebral regions as well as radiography of the extremities.

### 3.5.3 Operating elements


Operating elements can be found on both sides.

- (1) Stand column with LED indicator (optional) on the top of the column  
LED green: System is ready for exposure.  
LED yellow: Exposures are being performed.  
LED blue: System is not ready for exposure.
- (2) Emergency STOP push button
- (3) Front panel with imprint for detector position and measuring fields
- (4) Patient hand grips
- (5) Buttons for motorized vertical movement
- (6) Button for parking position
- (7) Button for centering tube to detector
- (8) Button for tube/detector tracking function with LED indicator
- (9) Button for manual vertical movement
- (10) Button for unlocking the tilt mechanism
- (11) Degree of tilting angle

### 3.5.4 Front panel markings


Example markings on wall stand with MAX wi-D (left) and Core XL / Core static (right)

- (1) Center crosshairs
- (2) Position of IONTOMAT ionization chambers
- (3) Detector field for MAX wi-D in portrait format
- (4) Detector field for MAX wi-D in landscape format
- (5) Center marking for portrait format
- (6) Detector field for Core XL / Core static

## 3.6 Wireless detectors (MAX wi-D or Core XL)

### 3.6.1 Description

- (1) Detector MAX wi-D
- (2) Detector Core XL
- (3) Detector Indicators
- (4) Wireless antenna

### General Description

The wireless detector is part of a digital image acquisition chain in an overall radiological system. It features a portable equipment designed for mobile applications.

The communication is managed through a WiFi interface.


Do not cover the wireless antennae during exposures, as otherwise the detector will lose its wireless signal.

## Detector description

The MAX wi-D and Core XL both have LED indicators for internal status and connector for service issues. The back sides (with respect to the active sensitive array) of the both detectors include electrical contacts to recharge the replaceable battery.

## Detector indicators (MAX wi-D)

- (1) Detector status LED
- (2) Wi-Fi status LED
- (3) Battery status LED
- (4) Infrared sensor  
needed for attachment of the detector
- (5) Power button  
Long push (switch off/on)

It may be possible that even when the Wi-Fi Status LED and the detector Status LED are green, release of radiation is blocked by the system.


Please check information displayed on the user interface (TUI).

LEDs are provided on the detector (close to handling), to give basic information on the detector behavior.

| LED color / status         | Detector status LED (1)                             | Wi-Fi status LED (2)                                        | Battery status LED (3) |
|----------------------------|-----------------------------------------------------|-------------------------------------------------------------|------------------------|
| LED off                    | OFF state                                           | Wi-Fi disabled                                              | No battery             |
| LED green                  | Operating state and online<br>Ready for acquisition | Wi-Fi available (detector connected to an Access Point)     | Battery OK             |
| LED orange                 | Offline or download state                           | Wi-Fi not ready (detector not connected to an Access Point) | Battery getting low    |
| LED orange - blinking slow | Error state or switching ON                         | n.a.                                                        | n.a.                   |
| LED orange - blinking fast | n.a.                                                | n.a.                                                        | Battery low            |

| LED color / status                                                      | Detector status LED (1)                                                     | Wi-Fi status LED (2) | Battery status LED (3) |
|-------------------------------------------------------------------------|-----------------------------------------------------------------------------|----------------------|------------------------|
| LED orange - blinking fast, then red flashlight on the 3 LEDs, then OFF | Switching OFF (going to OFF state) when requested by message or by time-out | n.a.                 | n.a.                   |

### Detector indicators (Core XL)

- (1) Detector status LED
- (2) Wi-Fi status LED
- (3) Battery status LED
- (4) Cable port  
needed for attachment of the detector only

It may be possible that even when the Wi-Fi Status LED and the detector Status LED are green, release of radiation is blocked by the system.

Please check information displayed on the user interface (TUI).

LEDs are provided on the detector (close to handling), to give basic information on the detector behavior.

| LED color / status          | Detector status LED (1)          | Wi-Fi status LED (2)                                        | Battery status LED (3)       |
|-----------------------------|----------------------------------|-------------------------------------------------------------|------------------------------|
| LED off                     | n.a.                             | n.a.                                                        | No battery                   |
| LED green                   | Detector connected to system     | Wi-Fi available (detector connected to an Access Point)     | Battery capacity: 11% - 100% |
| LED orange                  | Detector not connected to system | Wi-Fi not ready (detector not connected to an Access Point) | Battery capacity: 6% - 10%   |
| LED orange - blinking       | Detector in firmware update      | n.a.                                                        | Battery capacity: ≤ 5%       |
| LED green - blinking        | Detector in startup              | n.a.                                                        | Battery charging             |
| LED green - blinking (once) | System in exposure               | n.a.                                                        | n.a.                         |

### 3.6.2 Status and system messages

In the following situations, the battery will probably be empty or the voltage will be too low:

- The wireless detector has been accidentally left outside one of the charging locations for too long, for example over night.
- The wireless detector is outside the room.

Proceed as indicated in the system messages to solve the problem.

## 4 System Operation

### 4.1 On-Off

The MULTIX Impact C system may be in one of the following operating status:

- Off

All the components of the system and the imaging station are turned off except for the detectors which remain in standby mode.

- On

All the components of the system are turned on and may be supplied from the power line. When the start-up phase is completed, examinations and post-processings may be performed.

- Emergency SHUTDOWN

All the components including the detector power supply and the imaging system are disconnected from the mains using the on-site emergency SHUTDOWN switch.

#### 4.1.1 Before using the system

If the power network at the site is very unstable, it may be advisable to install an uninterruptible power supply to prevent data loss. Contact Siemens Healthineers for more information.

- 1 Perform a visual inspection on the system.
- 2 Make sure that there is nothing to hinder stand movements.

#### 4.1.2 Switching on a switched-off system

No indicators on the generator ON/OFF console are lit.

- ◆ Press the ON button.

The ON indicator is lit.

The system is powered up.

The system performs a self-test.

The imaging system is started.

If the battery of the detector is completely empty, charge it and wait at least 60 minutes. The detector needs a warm-up time between 60 and 120 minutes to achieve good image quality, depending on environmental temperature and switch-off time.

During warm-up time acquisitions are possible but reduced image quality may occur.

(This is the case if the system has been turned off completely via the mains switch and has not been in standby mode. If it has been turned off at the On/Off control console, no warm-up time is necessary.)

In case of a fault:

- ◆ Shut down the system and turn it on again.

If the fault remains:

- ◆ Call Siemens Healthineers Customer Service.

### 4.1.3 Switching on after power failure or after an emergency SHUTDOWN switch is activated

When the on-site emergency SHUTDOWN switch is activated, the whole system is cut-off from the power supply.

Unstored images of the patient currently being examined may be lost. Images and data of terminated examinations of previous patients are stored to the hard disk of the image processor.

- 1 After a power failure or an emergency SHUTDOWN switch is activated, check the images of the last examined patient!
- 2 After removing the reason for activating the emergency SHUTDOWN switch, deactivate the emergency SHUTDOWN switch and then switch on the system.
- 3 If the imaging system does not start any more or if a serious error is indicated, call Siemens Healthineers Customer Service.

### 4.1.4 System switch-off

If no more acquisitions are required, you can switch the system off completely.

Do not press the keys **Ctrl+Alt+Del**.

Pressing this key combination is on your own risk.

If you do so, problems may occur.

- 1 Press the OFF button on the console.

The MULTIX Impact C system and the imaging system automatically shut down.

If the ON button is pressed again within 10 seconds (configurable by Siemens Healthineers Customer Service), only the MULTIX Impact C system shuts down, not the imaging system.

- 2 Switch off the main switch of the x-ray unit (or room), if necessary.

### 4.1.5 Emergency SHUTDOWN switch activation

- ◆ Only (exclusively) if patients, operators, third parties, or the unit are in danger, press (activate) the on-site emergency SHUTDOWN switch.

### 4.1.6 Rebooting the system

We recommend a daily reboot of the system or at least one reboot a week.

- 1 Press the OFF button on the console.
- 2 When the Image Station has shut down (dark monitor), press the ON button again.

### 4.1.7 UPS (for the imaging system)

UPS = Uninterruptible Power Supply

A UPS is a battery-backed system which provides emergency power in the event of a mains power supply failure.

To prevent data/image loss, your system may have a UPS which shuts down the imaging system in a controlled way in the event of a power failure.

#### CAUTION

Inappropriate device connected to UPS

**Risk of electrical hazard or damage to the system**

- ◆ Do not connect any additional device to the UPS power outlets.


### Functioning in the event of power failure

In the event of power failure, all indicators and displays are dark. Only the imaging system as well as the main monitor in the control room is still being powered by the UPS to end all active tasks, for example archiving of images.

- 1 Stop the examination.
- 2 Do not press any keys.  
After a configurable time (typically 2 minutes), the imaging system shuts down to prevent the backup battery to be exhausted.
- 3 Wait for power to return or start emergency power.
- 4 Restart the system.


### Description of UPS

The UPS has a five-button graphical LCD. It provides useful information about the UPS itself, load status, events, measurements and settings.

- (1) Power ON indicator (green)  
The UPS is operating normally
- (2) Battery ON indicator (yellow)  
The UPS is on battery mode
- (3) Alarm indicator (red)  
The UPS has an active alarm or fault
- (4) Escape
- (5) Up
- (6) Down
- (7) Enter
- (8) ON/OFF button

#### LCD description

- (1) Operation status
- (2) Load/equipment status
- (3) Efficiency and load group information
- (4) Battery status

As default, or after 5 minutes of inactivity, the LCD displays the screen saver.

The backlight LCD automatically dims after 10 minutes of inactivity.

Press any button to restore the screen.

The following table describes the status information provided by the UPS:

| Operation status                                                                                                                                 | Possible cause                                                                           | Action                                                                                                                                                                                             |
|--------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Standby mode<br>                                                | The UPS is OFF, waiting for start-up command from user                                   | System is not powered until power button is pressed.                                                                                                                                               |
| Normal mode<br>                                                 | The UPS is operating normally.                                                           | The UPS is powering and protecting the system.                                                                                                                                                     |
| In AVR mode<br><br>No beep                                      | The UPS is operating normally, but the system voltage is outside normal mode thresholds. | The UPS is powering the system through an Automatic Voltage Regulation device. The system is still protected normally.                                                                             |
| On battery<br><br>Battery LED is on<br>1 beep every 10 seconds | A utility failure has occurred and the UPS is in battery mode.                           | The UPS is powering the system with the battery power. Prepare your system for shutdown.                                                                                                           |
| End of backup time<br><br>1 beep every 3 seconds              | The UPS is in battery mode and the battery is running low.                               | This warning is approximate and the actual time to shutdown may vary significantly.<br><br>Depending on the UPS load, the "Battery low" warning may occur before the battery reaches 20% capacity. |

### Operation

#### Start-up and normal operation

To start the UPS:

- 1 Verify that the UPS power cord is plugged in.  
The UPS front panel display illuminates and shows EATON logo.
- 2 Verify that the UPS status screen shows .
- 3 Press the  button on the UPS front panel for at least 2 seconds.  
The UPS front panel display changes status to "UPS starting...".
- 4 Check the UPS front panel display for active alarms or notices. Resolve any active alarms before continuing.  
(→ Page 87 Troubleshooting)

If the indicator is on, do not proceed until all alarms are clear.  
Check the UPS status from the front panel to view the active alarms.  
Correct the alarms and restart if necessary.

- 5 Verify that the indicator illuminates solid, indicating that the UPS is operating normally and any loads are powered and protected.  
The UPS should be in normal mode.

### Starting the UPS on battery

Before using this feature, the UPS must have been powered by utility power with output enabled at least once.

To start the UPS on battery:

- 1 Press the button on the UPS front panel until the UPS front panel display illuminates and shows the status "UPS starting...".  
The UPS cycles through Standby mode to Battery mode.  
The indicator illuminates solid.  
The UPS supplies power to your system.
- 2 Check the UPS front panel display for active alarms or notices besides the "Battery mode" notice and notices that indicate missing utility power.
- 3 Resolve any active alarms before continuing.  
See (→ Page 87 Troubleshooting)
- 4 Check the UPS status from the front panel to view the active alarms.
- 5 Correct the alarms and restart if necessary.

### UPS shutdown

To shut down the UPS:

- ◆ Press the button on the UPS front panel for 3 seconds.  
The UPS starts beeping and shows the status "UPS shutting OFF...".  
The UPS then transfers to Standby mode and the indicator turns off.

### Operation on battery power

Transfer to battery power:

- The connected devices continue to be supplied by the UPS when AC input power is no longer available.  
The necessary energy is provided by the battery.
- The and indicators illuminate solid.
- The audio alarm beeps every 10 seconds.

The connected devices are supplied by the battery.

Low-battery warning:

- The  and  indicators illuminate solid.
- The audio alarm beeps every 3 seconds.

The remaining battery power is low.

Shut down all applications on the connected system because automatic UPS shutdown is imminent.

End of battery backup time:

- LCD displays "End of backup time".
- All LEDs go OFF.
- The audio alarm stops.

Return of AC input power:

- Following an outage, the UPS restarts automatically when AC input power returns and the load is supplied again.

### Troubleshooting

| Operation status                                                                                              | Possible cause                                                                                                     | Action                                                                                                                                                                                                                              |
|---------------------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Batteries disconnected<br> | The UPS does not recognize the internal batteries                                                                  | If the problem persists, contact Siemens Healthineers service.                                                                                                                                                                      |
|                                                                                                               | The batteries are disconnected                                                                                     | Verify that all batteries are properly connected.<br>If the problem persists, contact Siemens Healthineers service                                                                                                                  |
| Overload<br>               | Power requirements exceed the UPS capacity (greater than 105% of nominal)                                          | Remove some of the devices from the UPS. <ul style="list-style-type: none"> <li>• The UPS continues operation, but may shutdown if the load increases.</li> <li>• The alarm resets when the condition changes to normal.</li> </ul> |
| End of battery life<br>    | The end of the battery life is reached.                                                                            | Contact Siemens Healthineers service for battery replacement.                                                                                                                                                                       |
| Event<br>                  | A UPS event occurs.                                                                                                |                                                                                                                                                                                                                                     |
|                                                                                                               | Example:<br><br>Remote power OFF, the RPO contact has been activated to shutdown the UPS and now prevents restart. | Set the contact back to its normal position and press  button to restart.                                                                      |

| Operation status                                                                               | Possible cause                 | Action                                                                                                                                           |
|------------------------------------------------------------------------------------------------|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------|
| UPS fault<br> | The UPS has an internal fault. | The UPS does not protect the system anymore.<br><br>Note the alarm message and the UPS serial number, then contact Siemens Healthineers service. |

## 4.2 Functional and safety check

### 4.2.1 Daily checks

#### After system switch-on

- ◆ Perform a visual inspection of all displays and signal lamps on the operating units.  
No errors must be indicated.  
If radiation is released, all radiation ON indicators which may be present must light up.

#### Damage check

- 1 Check system and accessories daily for damaged or loose parts.
- 2 If problems are found, have the parts repaired before resuming patient examinations.

#### Prior to examination

- 1 Remove unnecessary objects and equipment from the area of action of the system.
- 2 Remove unnecessary accessories from the accessory rails of the table and the primary collimator.
- 3 Attach the devices required for patient positioning and immobilization securely to the system.
- 4 Attach all safety-related accessory parts correctly (for example, footboard, grip protection strip, hand grip, hand grip strip) and make sure that they are properly secured.
- 5 Clean any contrast agent residues from the patient table, protective plate.

#### Recommendation:

- 1 As a test, take a full-format exposure using the largest possible format.
- 2 Run a functional test of the emergency STOP buttons by performing an arbitrary system movement and pressing the emergency STOP button during the movement.  
The movement concerned must stop immediately.
- 3 Then unlock the button turning it clockwise.  
All system movements are possible again.

Check functionality of emergency stop daily.

### Checking system movements

#### Overhead support

See also (→ Page 94 *Positioning the 3D overhead support with tube unit*).


#### WARNING

Abnormal sinking of the overhead column.

#### Injury of persons

- ◆ Perform maintenance at regular intervals (at least every 24 months) to maintain the safety and proper functioning of the product.


- 1 Check manual longitudinal, transverse and vertical movements of the overhead support.
- 2 Check 'floating' overhead support movement.
- 3 Check the unit rotation.
- 4 Check motor driven overhead support height adjustment.
- 5 Observe the indications in the display.

#### Patient table (optional)

- 1 Check manual longitudinal and transverse movement of the tabletop.
- 2 Check motor driven table height adjustment.

#### Wall stand

Verify up and down movement of the detector unit at the wall stand:

- ◆ Press the detector unit's up and down buttons.

#### Checking the light localizer and centering/collimation

Check the center position of radiation unit and wall stand with the light localizer and check the collimation.

#### Light on

- ◆ Press this button at the multileaf collimator.  
The light beam for collimation appears (time limited).

### Centering

The accuracy of the stop positions of the tube unit can decrease to  $\pm 3$  degrees due to hard wear. Lateral image information may be cut off.

Always check the alignment of the tube unit to the detector with the light localizer.



- ◆ Check the alignment of the tube unit to the detector with the light localizer.

### Collimation

Collimation close to the object reduces the scattered radiation and thus improves image quality. The radiation exposure to the patient is also reduced.

To protect the patient from unnecessary radiation, always check the collimator position with the light localizer.

- ◆ Open and close diaphragm.  
The light beam changes its height or width.

Checking the automatic collimation (ACSS): When ACSS is displayed on the TUI, the automatic collimation is active.

### Light off

- ◆ Press the button again if light is on.  
The light beam is turned off.

### During examination

- 1 Check the radiation ON indicator.

It should light up or sound only if one of the radiation release switches has been pressed or for the duration of an X-ray exposure.

#### **CAUTION**

X-ray does not stop due to a technical defect or unintentional exposure. The radiation indicator or a radiation warning lamp lights up without a button/switch being activated.

#### **Risk of unnecessary radiation exposure**

- ◆ Press the SHUTDOWN button to switch off the whole system.
- ◆ If no cause can be found, contact Siemens Healthineers Service and discontinue use of the system.

- 2 Check the patient positioning devices, for example the grip protection strip and hand grips.
- 3 Initiate system movements only under the following conditions:
  - There is no danger of injury to the patient or third parties.
  - There are no objects blocking the path of system movements.

### 4.2.2 Monthly checks

Make the following checks in addition to the daily checks:

#### SID

The three-digit digital SID display = source-image distance can be checked for agreement with the tape measure in the multileaf collimator.

We recommend you to make this check at regular intervals - approximately monthly.

Incorrect SID values can cause overframing.

In this case please call the Siemens Healthineers Service.

### Constancy test

Perform the constancy test according to the applicable (country-specific) national statutory provisions.

#### CAUTION

Image incorrectly exposed because of malfunction of AEC

#### Unexpected radiation

- ◆ Perform a functional test of the automatic exposure control system as described below.

### 4.2.3 Functional check of the automatic systems

#### Automatic exposure control

If the values obtained in the following tests deviate from those specified, switch the system off and contact Siemens Healthineers Service immediately.

- 1 Select the an organ program.
- 2 Close the collimator on the tube assembly to a format of 10 cm x 10 cm.
- 3 Place a lead rubber apron, folded four times, in the beam path as a cover.
- 4 Release an exposure and keep the exposure switch pressed down.

The automatic exposure control automatically switches off when the limit is reached.  
A message appears: "Acquisition abort on limit excess".

- 5 Release the exposure switch.
- 6 Open the collimator on the tube assembly and remove the lead rubber apron.
- 7 Press the exposure switch.

The automatic exposure control switches off automatically at  $< 1$  mAs.

#### 4.2.4 Checks concerning radiology

- IQAP check at delivery
- Maintenance every 24 months
- Daily image check by customer, for example, for artefacts.


#### 4.2.5 DAP display test

The phantom for the test of the dose area product display consists of a 20 cm thick rectangular block of polymethylmethacrylate (PMMA) with side lengths greater than or equal to 25 cm.

- 1 Measure the air kerma at a distance of 60 cm from the focal spot (approx. 55 cm above the tabletop for an SID of 115 cm) and multiply this dose value by the cross-sectional area of the useful beam in this plane.
- 2 Convert the result into units of  $[\mu\text{Gy} \times \text{m}^2]$  and compare it with the value displayed.
- 3 In case of inconsistent DAP values, call service.

### 4.3 Preparation and movements of the overhead support

#### 4.3.1 Directions for moving and rotating the tube

Example

#### 4.3.2 Directions for moving

The 3D ceiling-mounted tube support can be moved in the following directions:

- x direction: longitudinal movement
- y direction: transverse movement
- z direction: height movement

#### 4.3.3 Directions for rotating

The X-ray tube unit can be rotated in the following directions:

- about the horizontal axis (RHA)
- about the vertical axis (RVA)

#### 4.3.4 Rotating the collimator

For more detailed information refer to (→ Page 99 *Rotating the multileaf collimator*)

### 4.3.5 Positioning the 3D overhead support with tube unit

#### WARNING

Defects on the carrying ropes or spring inside the lifting column broken

##### **Injury of persons**

- ◆ Do not exert strong force.
- ◆ Call the Siemens Healthineers service.

#### WARNING

Damaged cables inside the lifting column of the overhead tube or detector.

##### **Overhead tube or detector can fall down. Risk of patient or user injury.**

- ◆ If the overhead tube or detector cannot be moved vertically, do not force the movement.
- ◆ Call Siemens Healthineers Service.

#### CAUTION

In extreme operating mode, the X-ray tube assembly can heat up.

##### **Danger of burns**

- ◆ Avoid contact with the tube housing to avoid burns.

Detents are provided in the longitudinal and transverse direction of the 3D overhead support, i.e. the support can be moved only by applying increased force or by releasing the button and pressing it again to go through certain positions (e.g. table center in longitudinal and transverse direction or defined SID values in longitudinal direction).

The buttons in the handles are configurable by the Siemens Healthineers Uptime Service:

- Floating movement in X-, Y- and Z-direction
- Floating movement in X-, and Y-direction
- Floating movement in Z-direction
- With or without detent functioning in X- and Y-direction

### Adjusting height (z direction) manually

- 1 Grasp the handles of the tube unit with both hands
- 2 Press the button to release the brake for vertical movement.  
Releasing the button will stop the movement.
- 3 Pull the tube unit into the required position.

### Adjusting height (z direction) motor driven

- ◆ Press one of the buttons on the control panel for up or down movement.

### Adjusting longitudinal position (x direction)

- 1 Grasp the handles of the tube unit with both hands
- 2 Press the button to release the brake for longitudinal movement.  
Releasing the button will stop the movement.
- 3 Pull the tube unit into the required position.

Detents are provided in the longitudinal and transverse direction of the overhead support, i.e. the support can be moved only by applying increased force or by releasing the button and pressing it again to go through certain positions (e.g. table center in longitudinal and transverse direction or defined SID values in longitudinal direction: 115 cm, 150 cm, 180 cm and 300 cm for Ortho acquisition).

### Adjusting transverse position (y direction)

- 1 Grasp the handles of the tube unit with both hands
- 2 Press the button to release the brake for transverse movement.  
Releasing the button will stop the movement.
- 3 Pull the tube unit into the required position.

Detents are provided in the longitudinal and transverse direction of the overhead support, i.e. the support can be moved only by applying increased force or by releasing the button and pressing it again to go through certain positions (e.g. table center in longitudinal and transverse direction or defined SID values in longitudinal direction: 115 cm, 150 cm, 180 cm and 300 cm for Ortho acquisition).

### Rotating the tube unit (RHA/RVA)

The tube unit (with multileaf collimator) can be rotated manually about two axes:

- about the horizontal axis (RHA)
- about the vertical axis (RVA)

#### CAUTION

In extreme operating mode, the X-ray tube assembly can heat up.

#### **Danger of burns**

- ◆ Avoid contact with the tube housing to avoid burns.

The temperature of the tube surface during normal use reaches 50°C.  
In case of high energy radiation the temperature may exceed the 50°C.

#### Rotation about horizontal axis (RHA)

- 1 Grasp the handles of the 3D overhead support with both hands.
- 2 Press the button to release the brake for rotating the tube unit about the horizontal (support arm) axis with your right thumb.

- 3 Turn the tube unit into the required position.  
The angle of rotation is displayed in degrees on the **Position** tab card on the TUI.

The tube unit can be rotated about the horizontal axis continuously, at a maximum by  $\pm 140^\circ$ . There are stop positions at  $0^\circ$  and  $\pm 90^\circ$ .

When the software stop function is active:

- 1 Rotate the tube unit to a position  $\pm 3^\circ$  from the required detent ( $0^\circ$  or  $\pm 90^\circ$ ).
- 2 Remove your finger from the button but hold the tube handles for one second.

The tube unit will be rotated automatically to the detent.

The automatic rotation will be interrupted if any other tube movement happens in the meantime.

Press the RHA button again to reactivate the function.

(→ Page 64 *Software stop (only available for systems with Ortho or SmartPositioning function)*)

The accuracy of the stop positions of the tube unit might decrease to  $\pm 3^\circ$ , due to hard wear for example. Lateral image information might be "cut".

Always check the alignment of the tube unit with the light localizer.

#### CAUTION

Rotation of tube

##### Risk of crushing

- ◆ Pay special attention to all crushing hazards between the system's moving parts and their guide openings.

#### Rotation about vertical axis (RVA)

- 1 Grasp the handles of the 3D overhead support with both hands.


- 2 Press the button to release the brake for rotating the tube unit about the vertical (column stand) axis with your left thumb.

You can remove your thumb from the button once the brake has been released.

- 3 Turn the tube unit into the required position.

The tube unit can be rotated about the vertical axis from  $+182^\circ$  to  $-154^\circ$  with stop positions at  $+180^\circ$ ,  $+90^\circ$ ,  $0^\circ$  and  $-90^\circ$  (see drawing).

Rotation of the tube unit about the vertical axis (view from above: from the ceiling to the floor)

#### 4.3.6 Tube tracking

##### Tube tracking in table acquisition mode

Prerequisite:

- Workplace "Patient table" is selected and a detector is inserted.

Press the button on the TUI to activate the tracking function.

Press the button for rotation about the horizontal axis on the tube and rotate the tube into the required position.

- or -

Press the button for longitudinal movement on the tube and pull the tube into the required position.

With RHA of the tube up to  $45^\circ$  in both directions, the detector inside the table follows the tube rotation or longitudinal movement; the central beam is locked on the detector.

The detector unit stops movement when it reaches the required position or the button for tube movement has been released.

##### Tube tracking in wall acquisition mode

Prerequisite:

- Workplace "Wall stand" is selected and a detector is inserted.

Press the button on the wall stand or TUI to activate the tracking function.  
The LED indicator on the wall stand is on.

Press one of the buttons for vertical movement on the wall stand to move the detector up and down to the required working height.

With RHA of the tube to  $90^\circ$ , the tube follows the detector inside the wall stand vertically; the central beam is locked on the detector.

The tube stops movement when it reaches the required height or the button for vertical movement on the wall stand has been released.

The tube movement stops immediately when it reaches the collision zone.

#### 4.3.7 SID tracking

In addition to the tracking control perpendicular to the beam direction (in x direction), a further tracking control (in z direction) can be activated for keeping the SID (source-image distance) constant.

#### CAUTION


Parts of the system move unexpectedly during automatic positioning or tracking

#### Risk of injury to the patient or user

- ◆ Keep the patient within visual contact whenever performing any positioning functions.
- ◆ Use the emergency stop if motion needs to be stopped quickly.

#### Activating SID tracking in table acquisition mode

- ✓ Workplace table mode is active.
- ✓ A detector is inserted in the closed tray of the patient table.
- ✓ Tube unit is positioned above the table.
- ✓ Tube unit RHA angle is  $0^\circ \pm 3^\circ$ .
- ✓ Tube unit RVA angle is  $0^\circ$ .
- ✓  $SID \geq 90$  cm
- ◆ Touch the SID tracking button at the TUI.

### Activating SID tracking in wall acquisition mode

- ✓ Workplace wall mode is active.
- ✓ A detector is inserted in the closed tray of the wall stand.
- ✓ Detector unit is in  $+90^\circ$ /horizontal position.
- ✓ Tube unit is positioned above the detector unit.
- ✓ Tube unit RHA angle is  $0^\circ \pm 3^\circ$ .
- ✓ Tube unit RVA angle is  $0^\circ$ .
- ✓  $SID \geq 90$  cm
- ◆ Touch the SID tracking button at the TUI.

### 4.3.8 Deactivation of a tracking control

If the state of the system is changed so that one or several of the conditions for active tracking control are not fulfilled, then the tracking control is automatically switched off.

### 4.3.9 SID display

The SID (source-image distance) is displayed (in cm or inch according to configuration) at the TUI.

During oblique projections the angle in degrees is displayed on the TUI.

You must therefore absolutely measure the SID with the tape measure in the multileaf collimator.

### 4.3.10 Rotating the multileaf collimator


#### CAUTION

Hands between the handles and the overhead column

#### Risk of injury

- ◆ Grasp the collimator on the sides when turning it so that there is no risk of crushing your hands.

The collimator can be rotated out of the stop position  $0^\circ$ . The maximum rotation in both direction is  $\pm 45^\circ$  by using the Siemens Healthineers tube flange. Assure that the rotation is not limited by the wiring harness.

#### Rotating the collimator

- (1) Lock-in lever (not available on the Collimator RFU)
- (2) Lamp housing

If the collimator is NOT in the 0° position, it will be indicated at the display of the TUI with this symbol.

#### Rotating the collimator $\pm 45^\circ$ around the Central Beam Axis

- 1 Press the lock-in lever on the multileaf collimator in the direction of the front panel.  
The 0° lock of the multileaf collimator is released.  
(There is no 0° lock in the Collimator RFU.)
- 2 Grasp the collimator and rotate it by the required angle to the required direction.

#### Rotating back into the 0° lock-in position

- ◆ Grasp the collimator and turn it into the centered 0° lock-in position.

## 4.4 Preparation and movements of the patient table

### WARNING

System movements cause collision

#### Risk of injury or system damage

- ◆ Perform system movements only when it is certain that neither the operator, the patient, third parties nor pieces of equipment can be endangered by these movements.
- ◆ Remove any objects or accessories, for example injector or infusion stand, from the collision area.
- ◆ Make sure you are standing outside the dangerous area.
- ◆ Always watch the patient when performing system movements.

### WARNING

Brakes on tabletop have reduced force when the system is turned off.

#### Risk of injury

- ◆ Check that nobody is leaning on the tabletop when the system is turned off.
- ◆ Avoid positioning patients such that a sudden movement of the tabletop could cause them to fall.
- ◆ Execute patient transfer with special care.

The table can be moved with the foot kick switches in the table base.

Optionally, there is also a foot switch for the elevating patient table.

### 4.4.1 Collision protection

The table has a collision detection that stops downward motorized movements in case the tabletop collides with an obstacle.

The tabletop can nevertheless be moved upwards.

### 4.4.2 Adjusting table height

### WARNING

Lowering the table onto a seated patient

#### Crushing of patient or operator seated with legs under the table (for example, a wheelchair patient)


- ◆ Do not position a patient with legs under the table.
- ◆ Do not move the table while a patient is seated near the table. System will stop when collision is detected, but there is still a risk of injury for the patient!
- ◆ Move the patient away from the table first.

### Moving the table up

- ◆ Actuate the second foot kick switch from the right in the table base.

– or –

Actuate the middle foot switch (optional) for table up movement.


### Moving the table down

- ◆ Actuate the second foot kick switch from the left in the table base.

– or –

Actuate the left foot switch (optional) for table down movement.

The table stops automatically at a certain height which can be programmed by Siemens Healthineers service.

#### 4.4.3 Positioning the tabletop

##### **CAUTION**

Patient does not use the hand grips and gets fingers caught in moving parts.


##### **Risk of crushing**

- ◆ Make sure that patient always uses the hand grips.
- ◆ If this is not possible, make sure that the patient's hands cannot get caught in the moving parts of the system.

- 1 Actuate one of the two outer foot kick switches in the table base.

– or –

Actuate the right foot switch to release the brakes for tabletop floating.

The tabletop brakes are released.


- 2 Push the tabletop manually into the required position.

Putting a strain on the extended tabletop may activate the collision control device.

#### 4.4.4 Inserting wireless detector

- 1 Push the tabletop away from you as far as possible.

##### Example

- 2 Press the brake handle in the handgrip (A) to release the detector tray (B).
- 3 Pull out the detector tray (B) until it reaches the stop position.

##### CAUTION

Incorrect insertion of detector into the Bucky tray.

**Risk of pinching or damage to detector.**

♦ Be careful when inserting detector in the Bucky tray.

- 4 Carefully place the rear edge of the detector onto the detector tray, push the rear edge of the detector against the stops and lower the detector onto the tray between the lateral stops.

Only the MAX wi-D can be placed either in portrait or landscape format.

- 5 Carefully and slowly slide the detector tray into the patient table.

Push the detector tray back completely, or else the detector is not loaded in the detector tray.

Be careful not to damage the detector when sliding it into the tray as the tray does not have a mechanism to absorb the shock.

#### CAUTION

Detector cover broken

##### Toxic contamination

- ◆ Handle portable detector with care!
- ◆ If the detector is broken do not use it.
- ◆ Inspect if any particles are spilled.
- ◆ Leaked particles could be toxic - do not touch them!
- ◆ Use protective gloves to collect and keep them in a sealed container.
- ◆ Return spilled particles to Siemens Healthineers.

Do not load the withdrawn detector tray with more weight than the label on the detector tray states.

#### CAUTION

Detector tray overloaded

##### Risk of damage to the system and injury of the patient

- ◆ The maximum load of the opened detector tray must not exceed 10kg.
- ◆ Do not allow the patient to use the detector tray for support.

Do not pull the detector tray out of the patient table immediately after the exposure, wait at least 3 s after the exposure.

### 4.4.5 Positioning the detector tray

The detector tray can be moved manually to the left or right. The detector tray also follows the tube unit motor-driven within a certain range (pull-in range).

#### Longitudinally

- For normal exposures on the patient table, the detector tray can be moved longitudinally to the left or right using the break knob. To unlock the brakes you have to pull the knob.
- When the tracking function is active, the detector tray follows the longitudinal movement of the X-ray tube.

## Centering

- ✓ If the software stop function is configured, activate the function and rotate the tube unit to the software detent at 0°. (→ Page 96 *Rotation about horizontal axis (RHA)*)

- ◆ To center the detector (tray) to the tube unit press the centering button on the tube unit.

The detector unit stops movement when it reaches the center position or the centering button has been released.

## Free exposures

(applicable only for table with mobile detector)

For exposures without grid and without automatic exposure control, and for better accessibility, for example during exposures of the extremities, the detector tray can be pulled out like a drawer.

### Pulling out the detector tray

- ◆ Pull the detector tray out to the front up to the stop.

The detector tray is locked in the end position.

Gently let slide the detector tray into the end positions, because strong vibrations may lead to impact on image quality.

### CAUTION

Detector tray overloaded

#### Risk of damage to the system and injury of the patient

- ◆ The maximum load of the opened detector tray must not exceed 10kg.
- ◆ Do not allow the patient to use the detector tray for support.

### Pushing the detector tray back in

- ◆ Push the detector tray back in until it locks into position.

Push the detector tray back completely, or else the detector is not loaded in the detector tray.

## 4.4.6 Scattered radiation grid

If necessary for pediatric exposures, the grid can be removed from or inserted into the beam path.

The responsibility for selecting and inserting the scattered radiation grid to be used for an exposure lies with the examiner.

The intended SID is color-coded on each table grid using the following color scheme:

| Colors for grid handle |                                                                                   |           | Grid parameters |          |
|------------------------|-----------------------------------------------------------------------------------|-----------|-----------------|----------|
| SID (cm)               | Grid                                                                              |           | Ratio           | Lines/cm |
| 115                    |  | dark blue | 13              | 92 & 40  |

Only the same type of grid can be shared from other system.

#### CAUTION

Grid is dropped or not handled properly

**Risk of invisible damage and impaired image quality**

- ◆ Handle the grid with special care.

#### CAUTION

Changing the grid

**Crushing between the housings of the detector case.**


- ◆ Pay attention to danger zones / danger points of crushing between the housing of the detector case.

#### CAUTION

Image made with incorrect grid status

**Poor image quality leading to repeated X-ray exposure**


- ◆ The expected grid status (inserted, removed) is indicated on the imaging system.
- ◆ Manually remove or insert the grid from or into the detector tray as required before making the X-ray exposure.

Top of grid (facing the tube unit)

### Removing the grid

- 1 Push the tabletop away from you.

- 2 Press the locking lever at the left side of the grip and, keeping it straight, pull the grid out towards the front.
- 3 Remove the grid.

## Inserting the grid

Please pay attention to the label which indicates the correct focus side and attach the grid correctly.

- 1 Hold the grid horizontally, labels facing up.

Please take care that the SID range of the grid is compatible to the selected SID.

- 2 Insert the grid into the accessory rails.
- 3 Holding it horizontally and straight, push the grid all the way in.  
When the grid has reached the end position, the locking lever engages.

## 4.5 Preparation and movements of the wall stand

### 4.5.1 Safety information

#### Information about holding facilities

#### CAUTION

Patient exerts too much weight on the patient hand grip on the wall stand.

#### Risk of injury to the patient and damage to the wall stand

- ◆ Do not let the patient hang on the patient hand grip.
- ◆ Do not exert more than 25kg weight on the patient hand grips.

#### CAUTION

Detector unit overloaded because patient is sitting on it or falling/leaning against it.

#### Risk of damage to the detector and injury of the patient

- ◆ Make sure that the detector is not mechanically overloaded as for example by patient sitting on it.

#### CAUTION

Loose patient hand grips on the wall stand

#### Mechanical or personal injury

- ◆ Check that patient hand grips are in right position when moving detector up or down.

**⚠ WARNING**

Parts of the system move unexpectedly during automatic positioning or tracking

**Risk of injury to the patient or user**

- ◆ Keep the patient within visual contact whenever performing any positioning functions.
- ◆ Use the emergency stop if motion needs to be stopped quickly.

**⚠ CAUTION**

In BWS, no grid or lighter grid is inserted, and button for manual vertical movement is pressed.

**Risk of injury to the patient or user due to unintentional movement**

- ◆ Hold the patient hand grip nearby when press the button for manual vertical movement.

**Information about system movements****⚠ WARNING**

Broken drive chain in wall stand

**Risks of injury when moving front unit up or down**

- ◆ Do not activate vertical movements when a broken drive chain is detected.

**⚠ WARNING**

Break of the chain

**High movement forces in vertical direction**

- ◆ Stop working with the component and call service.

**4.5.2 Detector unit movements**

The detector unit can be moved vertically using the respective buttons at the detector housing.

**Vertical movement (height adjustment) of the detector unit**

- ◆ Press the button for the manual vertical movement and move the detector unit up or down to the required working height.

When moving down the detector unit, be careful not to hit the ground.

– or –

Press one of the buttons for the motorized vertical movement to move the detector up and down to the required working height.

Downward movement stops at 12 cm distance to floor.

You can continue movement manually.

### Tilting the detector unit

The detector unit can be continuously tilted to the required position between 0°/vertical and +90°/horizontal or up to -20°.

- ◆ Press the button for unlocking the tilt mechanism and tilt the detector tray to the required position.

#### CAUTION

Detector unit overloaded because patient is sitting on it or falling/leaning against it.

#### Risk of damage to the detector and injury of the patient


- ◆ Make sure that the detector is not mechanically overloaded as for example by patient sitting on it.

### 4.5.3 Inserting the detector into the wall stand

(with wireless mobile detector only)

- 1 Move the detector unit to a convenient working height.


Vertical movement (height adjustment) of the detector unit

Example, detector tray for Core XL (left) and MAX wi-D (right)


- 2 Press the brake handle (A) and pull out the detector tray (B) until it reaches the stop position.

Example, detector tray with Core XL

- 3 Insert the detector from above into the holders as shown in above figure (with the centering cross facing toward the front).

Example, detector tray with MAX wi-D

Only the MAX wi-D can be inserted either in portrait or landscape format. When you insert the MAX wi-D in landscape format (see above figure on the left) you might have the feeling that the detector is not fixed properly. But be assured that the detector is fixed by the clamp (C) top right.

- 4 Press the brake handle.
- 5 Carefully and slowly slide in the detector tray with the detector.

Be careful not to damage the detector when sliding it into the tray as the tray does not have a mechanism to absorb the shock.

#### 4.5.4 Removing the detector

(with wireless mobile detector only)

- 1 Move the detector unit to a convenient working height.  
Vertical movement (height adjustment) of the detector unit
- 2 Use the brake handle (A) to pull out the detector tray with the detector.
- 3 Remove the detector pulling it upwards.
- 4 Carefully and slowly insert the detector tray again.

### 4.5.5 Scattered radiation grid

The responsibility for selecting and inserting the scattered radiation grid to be used for an exposure lies with the examiner.

The intended SID is color-coded on each wall stand grid using the following color scheme:

| Colors for grid handle |                                                                                    |               | Grid parameters |          |
|------------------------|------------------------------------------------------------------------------------|---------------|-----------------|----------|
| SID (cm)               | Grid                                                                               |               | Ratio           | Lines/cm |
| 115                    |   | dark blue     | 13              | 92 & 40  |
| 140*                   |   | grey          | 13              | 92 & 40  |
| 180                    |   | dark orange   | 13              | 92 & 40  |
| 300                    |  | bright yellow | 13              | 92       |

\* Universal grid (140), working SID from  $f_0=115$  cm to  $f_0=180$  cm

Only the same type of grid can be shared from other system.

#### CAUTION

Grid is dropped or not handled properly

**Risk of invisible damage and impaired image quality**

◆ Handle the grid with special care.

#### CAUTION

Changing the grid

**Crushing between the housings of the detector case.**

◆ Pay attention to danger zones / danger points of crushing between the housing of the detector case.


**CAUTION**

Image made with incorrect grid status

**Poor image quality leading to repeated X-ray exposure**

- ◆ The expected grid status (inserted, removed) is indicated on the imaging system.
- ◆ Manually remove or insert the grid from or into the detector tray as required before making the X-ray exposure.

### Removing the grid

- 1 Press the locking lever (A) and, keeping it straight, pull the grid out towards the side.
- 2 Remove the grid with both hands.

### Inserting the grid


When inserting the scattered radiation grid, make sure that the grid being inserted is correctly assigned to the intended source-image distance (SID) on which the grid is focused.

- 1 Hold the grid vertically, labels facing the tube unit.

The label with the tube symbol on the grid must point toward the tube unit.

- 2 Place the grid against the top and bottom guide rail.
- 3 Holding it vertically and straight, push the grid all the way in.  
When the grid has reached the end position, the locking lever engages.

#### 4.5.6 Positioning the tube unit automatically

- ✓ A suitable organ program for an exposure onto the wall stand is selected.
- ◆ Press the SmartPositioning button on the Wireless Remote Control Console.  
(SmartPositioning is an optional function.)  
The tube unit moves to the position provided in the selected organ program.

#### 4.5.7 Positioning the tube unit manually

- ◆ Position the tube unit manually as required for the exposure using the movement buttons at the TUI or the buttons on the handles.

When the tube unit has reached the correct longitudinal and transverse centering detents, the centering displays (white ring) at the movement buttons light up.

– or –

Activate the tube tracking function by pressing the button on the wall stand.

The LED indicator is on.

(→ Page 97 Tube tracking)

– or –

Move the tube unit to the vertical center position with the help of the light localizer.

### 4.5.8 Centering the tube unit in vertical direction

- ✓ If the software stop function is configured, activate the function and rotate the tube unit to the software detent at  $+90^\circ$  or  $-90^\circ$ . (→ Page 96 *Rotation about horizontal axis (RHA)*)

- ◆ Press the centering button on the wall stand.

The tube unit stops movement when it reaches the center position or the centering button has been released.

Centering is not possible when the tube unit reaches the collision zone.

– or –

Move the tube unit manually to the vertical center position with the help of the light localizer.

### 4.5.9 Moving the tube unit to the parking position

A parking position has been set as default for the user to easily reach the tube unit.

The parking position can be configured by Siemens Healthineers service.

- ◆ Press the button on the wall stand.

– or –

Press the Parking position button on the Wireless Remote Control Console.

## 4.6 Preparation and handling of the wireless detectors

### 4.6.1 General safety information

Do not install the detector close to facilities where water is used.

Do not spill liquid or chemicals onto the detector or, in cases where the patient is injured, allow it to become wet with blood or other body fluids, as doing so may result in fire or electric shock.

In such situation, protect the detector with a single use plastic bag.

Do not spray the detector directly for cleaning and disinfection purposes.

## Interference with other devices

The portable wireless detectors may be interfered with by other equipment, including portable and fixed RF communication equipment, even if such equipment meets the applicable emissions requirements.

## Image transfer may be interrupted sporadically.

- The operator must ensure that other wireless devices in the 2.4 GHz or 5 GHz ISM band should not be operated in the vicinity of the system.
- Please observe and verify normal operation of the portable wireless detectors prior to using it.

### 4.6.2 Recommendations

#### Temperature

##### CAUTION

If the system has been completely without power by mistake or due to a power failure, the detector will not be at its operating point.

##### **Poor detail resolution with all types of radiation!**

- ◆ Note that image quality is unreliable under such circumstances. It is the physician's responsibility to decide whether to continue work on the basis of such image quality.

If a detector temperature warning message appears, please check if the warm-up time from a cold start (120 minutes) is respected. If yes, we recommend a recalibration.

##### CAUTION

Detector temperature out of range

##### **X-ray without diagnostic value**

- ◆ Check the temperature icon on the display regularly.

## EMC compliance


Portable and mobile RF communication equipment can affect the detector operation.

## Warm-up time from a cold start

After a cold start (battery completely empty), it is possible to operate the detector, but with reduced image quality. Therefore, avoid the battery getting empty.


The detector reaches good image quality again after a maximum of 120 minutes (typically 60 minutes).

### 4.6.3 Notes on handling the detector



Handle the detector carefully, as it may be damaged if something is hit against it, if it is dropped from high altitude, or if it receives a strong jolt.

The detector is monitored by integrated sensors, so that data can be read out by a technician.

Be sure to use the detector unit on a flat place so it will not bend. Otherwise, the detector may be damaged.

Do not hit the detector surface in any way.

**⚠ CAUTION**

Detector cover broken

**Toxic contamination**

- ◆ Handle portable detector with care!
- ◆ If the detector is broken do not use it.
- ◆ Inspect if any particles are spilled.
- ◆ Leaked particles could be toxic - do not touch them!
- ◆ Use protective gloves to collect and keep them in a sealed container.
- ◆ Return spilled particles to Siemens Healthineers.

**⚠ CAUTION**

Detector cover breaks because patient is too heavy

**Possible injury of patient or personnel by carbon or glass splinters or electrical shock.**

- ◆ Use the FD protector for patients heavier than 100 kg (maximum 227 kg).
- ◆ Handle the portable detector with care!
- ◆ When the detector is broken, do not use it.

**⚠ CAUTION**

Hot surface of detector The surface of the portable detector can reach 43°C (109°F), which may be too hot for high risk patients, such as infants.

**Patient injury**

- ◆ Protect such patients from exposure to the surface of the detector.

**⚠ CAUTION**

Portable detector is too hot.

**Burns or heat injuries**

- ◆ Pay attention to the warning messages pertaining to detector surface temperature.
- ◆ Remove the detector from contact with the patient if the surface temperature exceeds 41°C (106°F).

**CAUTION**

Improper or careless detector handling

**Damage to the detector or detector surface**

- ◆ Do not drop the detector.
- ◆ Only use the detector on a level surface to avoid sagging.
- ◆ Exercise absolute care in avoiding pressure or sudden impact on the detector surface.

**CAUTION**


Exposure on a wireless detector that is connected to another system

**Image data will be lost.**

- ◆ Make sure that the detector and the system have the same colored labels.

**Shock sensor**

Both the MAX wi-D and Core XL have integrated shock sensors.

Core XL (back side)

The shock sensor (A) is located on the back side of the Core XL. Normally the shock sensor is white, and it turns red if the detector is dropped from high altitude (> 70 cm).

The shock sensor is not visible on the MAX wi-D. However, the shock details can be accessed in the **Detector Overview**.

(→ Page 126 FD management)

#### 4.6.4 Maintenance

No other maintenance action is required except the calibration needed for image quality, whose frequency is defined by the system program.

There are no parts that need to be replaced due to limited lifetime, except for the battery. In case of malfunction, the detector will be returned to the manufacturer for repair.

#### 4.6.5 Operation and setting of WLAN

##### WARNING

The portable detector generates an electromagnetic field so high that it disturbs life-supporting device.

##### Interference with life supporting device resulting in possible malfunction

- ◆ Keep a safety separation distance of more than 90 cm between the portable detector and the life-supporting device.
- ◆ In case of interference with other equipment increase the distance between the interfering devices.

##### CAUTION


Image cannot be made because the battery in the portable detector is low

##### Portable detector not available for use

- ◆ Connect the portable detector to a loading port when it is not in use.
- ◆ Consider purchasing the battery loading station and additional batteries if portable detector availability is important.

## Introduction

- (1) Base station (WLAN access point)
- (2) Wireless connection
- (3) Wireless detector
- (4) Imaging system
- (5) PACS

The wireless connection of the mobile detector is established between the X-ray detector and the system.

Therefore, the phase of wireless connection and transfer does not utilize the existing wireless network that may be installed in the clinic. However, some precautions should be taken that both can coexist without interferences.

## Properties of wireless connection

The system uses its own WLAN access point which is connected to the system and the X-ray detector as its sole client.

The communication is based on the WLAN standard 802.11. Within this standard, in most countries, there is the choice between the standards 11b, g, n (operating at 2.4 GHz) and the standards 11a, n (operation 5 GHz).

Within these standards (11a, b, g, n), several channels (frequencies) are available. A specific channel (frequency) can be selected by the site for use.

However, the available channels depend on legal regulations and requirements of the country with the installed system (wireless detector of MULTIX Impact C). The wireless connection is encrypted. For this the widely used and very secure WPA2 is applied.

The standard (11a, g, n) and the channel can be set by a Siemens Healthineers service technician via the service SW installed on the Image Station (imaging system). For a list of available channels in your country please contact your Siemens Healthineers project manager or sales representative.

The front end utilizes a power dependent on the standard and channel of the typical order of 32-79mW (always below 100mW). As the emission is not purely isotropic, more power is emitted in some directions and less in other directions. The maximum power emitted in a direction (max. EIRP) is 150 mW. This is defined based on the assumption that power would be emitted in all directions.

### Coexistence with other wireless connection

Wireless devices (e.g. your wireless hospital network and the mobile detector of the MULTIX Impact C system) can influence each other if they are operated at the same frequency.

For the mobile detector, this can lead to a slowdown or disruption of the wireless connection.

To avoid such a situation, it is necessary to enquire about wireless devices that are used in or within the vicinity of the installed MULTIX Impact C system. The contact person for enquiries regarding wireless network and devices inside your hospital regarding is probably your network or IT expert.

In case a WLAN-network is present in your work environment, the channels of the access points can be detected in the room where MULTIX Impact C will be installed. The frequencies in use (channels if you have WLAN) should be communicated to the Siemens Healthineers project manager.


Alternatively, select the channel to be utilized by the mobile detector of MULTIX Impact C and inform the service technician/project manager of the channel (frequency) accordingly.


It is important to ensure that this channel (frequency) will not be used at a later stage by other wireless devices.

### 4.6.6 Battery change (MAX wi-D)


The procedure is as follows:



- 1 Open the locking system of the battery as shown in the figure.
- 2 Remove the discharged battery.

- 3** Orientate the charged battery in the correct position: battery connector (1) facing the device connector (2).

- 4** Slide battery in the recess area of the device.  
**5** Push down the battery in the recess area.

- 6** Close the locking system of the battery and push down the latch to ensure no protrusion of the locking system.  
 When the battery is inserted correctly, the detector starts automatically.  
**7** Before using the detector, wait until all LEDs are green.

**CAUTION**


Incorrect handling of detector battery

**Explosion or fire**


- ◆ Handle the battery pack carefully.
- ◆ Do not drop, heat, open or mishandle it.

#### 4.6.7 Battery change (Core XL)


The procedure is as follows:

- 1 Open the locking system of the battery as shown in the above figure.
- 2 Remove the discharged battery.

- 3 Orientate the charged battery in the correct position: battery connector (1) facing the device connector (2).
- 4 Slide battery into the recess area of the device.
- 5 Push down the battery in the recess area.

- 6 Close the locking system of the battery.  
When the battery is inserted correctly, the detector starts automatically.
- 7 Before using the detector, wait until all LEDs are green.

#### **CAUTION**

Incorrect handling of detector battery

#### **Explosion or fire**

- ◆ Handle the battery pack carefully.
- ◆ Do not drop, heat, open or mishandle it.

### 4.6.8 Charging the battery of the detector

The battery of the detector can be charged either

- with the detector in the closed detector tray of the table\*, or
- with the detector in the closed detector tray of the Bucky wall stand\*, or
- in the battery charger

\*Automated detector battery charging as soon as the detector is located inside the table or Bucky wall stand tray, independent on the detector orientation.

Maximum operation time of a fully charged detector during regular utilization:

- MAX wi-D: 6.5 hours
- Core XL: 7.5 hours

Minimum charging time of a completely empty battery:

- MAX wi-D: 3 hours
- Core XL: 4 hours

Please put the detector in the closed detector tray for charging once acquisition has been finished.

When the system is switched off, make sure that the main power is on and put the detector in the closed detector tray, otherwise the battery will run down.


If the battery of the detector is completely empty, charge it and wait at least 60 minutes. The detector needs a warm-up time between 60 and 120 minutes to achieve good image quality, depending on environmental temperature and switch-off time.

#### 4.6.9 FD management

The wireless detector overview window at the imaging system shows the detector configuration of the system.

The wireless detector overview is opened by clicking **Detector Overview** in the control area of the **Patient** mode.

##### Example

- (1) It shows which type of detector is currently connected; the wireless detectors are shown with Wi-Fi, temperature and battery status
- (2) The serial numbers of the detectors help to identify the right detector.
- (3) It shows the date when the next calibration should be done.
- (4) It shows the current position of a detector.
- (5) It monitors the wireless detector and shows the shock history (only for the MAX wi-D).
- (6) The wireless detectors can be detached from the system for usage with another system.
- (7) Detector firmware version check  
If a new firmware version is available, click **Update** to install the new firmware. After the successful upgrade, the word **up-to-date** is displayed.

#### CAUTION

Wrong detector size connected or wrong detector used

##### Radiation not usable

- ◆ Check that detector identification on the detector label matches the detector identification shown on the user interface.
- ◆ Check that the green ready light on the correct detector is lighted before releasing radiation.

### 4.6.10 Attaching the detector

- ✓ The detector has been registered before to this system by the Siemens Healthineers service.
- ✓ Detector battery level > 50%
- ✓ **Detector Overview** dialog is displayed.

1 Connect the Core XL to the all-in-one PC with the cable.


– or –

Hold the MAX wi-D with the side of the LEDs to the IR sensor (1).

(1)

The detector can be taken away when the following message appears at the imaging system:

Example

- 2 Click **Yes** to attach the detector.
- 3 Wait a short time.

A message appears stating that the detector is attached successfully.

If the detector attachment was not successful, the detector has to be registered once to the system by the Siemens Healthineers Service.

#### Firmware upgrade during detector attachment

If a new detector firmware is available, the firmware upgrade will start automatically during detector attachment. Detector battery level < 50% will lead to firmware upgrade failure and the detector cannot be connected with the system.

- ◆ Charge the detector battery or replace it by a fully charged battery.
- ◆ It takes approximately 10 minutes to finish the firmware upgrade.
- ◆ Start the attachment again.

### 4.6.11 Detaching the detector

- ✓ The detector has been registered before to this system by the Siemens Healthineers Service.
- ✓ **Detector Overview** dialog is displayed.
- 1 Click **Disconnect** at the detector overview window.  
A message appears asking if you want to detach the detector.
- 2 Click **OK** to start the detachment procedure.  
A message appears stating that the detachment was successful.  
The detached detector is removed from the detector overview window.

### 4.6.12 Detector sharing with MAX wi-D

With the detector sharing function, it is possible to share the MAX wi-D detector between different systems (Multitom Rax, Ysio Max, Multix Fusion Max, Luminos dRF Max, Luminos Agile Max, MOBILETT Elara Max).

You can detach a MAX wi-D detector at the system where it is not needed and attach it to another system.

A detector calibration must be performed after the MAX wi-D has been exchanged.

Only one wireless detector can be connected to the MULTIX Impact C at the same time.

#### CAUTION

Exposure on a wireless detector that is connected to another system

**Image data will be lost.**

- ◆ Make sure that the detector and the system have the same colored labels.

### 4.6.13 Detector calibration

For details, please refer to the corresponding chapter in the imaging system operator manual.

**⚠ CAUTION**

Calibration procedure is started although there is somebody in the exam room.

**Unwanted radiation**

- ◆ Check if there is somebody in the exam room before starting the calibration.
- ◆ Do not leave the system unattended during calibration.
- ◆ If necessary, the calibration procedure can be easily stopped via the user interface.

**⚠ CAUTION**

Incorrect calibration parameters - Reduced image quality

**Risk of unnecessary X-ray exposure due to wrong calibration parameters**

- ◆ Make sure that the detector is recalibrated every 12 months to maintain the image quality.

**⚠ CAUTION**

Detector temperature out of range; Incorrect calibration data

**Risk of X-ray radiation without diagnostic value.**

- ◆ Check the temperature icon on the display regularly.

## 4.7 Emergency operation

For emergency operation (e.g., if the detector not working) you can also use the cassette for free exposures.

**⚠ CAUTION**

Portable detector cannot be connected to system due to hardware or software problems

**Portable detector not available for use**

- ◆ Maintain procedures and equipment so that cassettes can be used in place of the portable detector for free exposures if necessary.

## 5 Examination

### 5.1 Safety information for examination

#### WARNING

Brakes on tabletop have reduced force when the system is turned off.

##### **Risk of injury**

- ◆ Check that nobody is leaning on the tabletop when the system is turned off.
- ◆ Avoid positioning patients such that a sudden movement of the tabletop could cause them to fall.
- ◆ Execute patient transfer with special care.

#### WARNING

Hand grips loosen or are not used

##### **Injuries to patient from falling**

- ◆ Check that the patient hand grips on the table are in the right position.
- ◆ Ensure that hand grips are tightened prior to use.
- ◆ Make sure that the patient uses the hand grips during system movements or when turning over.
- ◆ Use other positioning aids for immobilizing weak or very unstable patients.

- 1 Attach all necessary safety accessories.
- 2 Use only positioning aids expressly approved by Siemens Healthineers. Contact your local Siemens Healthineers representative if you have any questions.  
  
Refer to (→ Page 157 *Accessories and Auxiliary Devices*) for all accessories available for positioning. Attaching the accessories is also described there.
- 3 Ensure that these devices and positioning aids are secured correctly and are functioning properly.
- 4 Make sure that the body parts of the patient, especially arms and legs, do not extend over the edge of the tabletop.
- 5 Make sure that the patient uses the provided grip locations.
- 6 Let the patient get on or off the table only in the central area.

#### CAUTION

Incorrect positioning of the patient.

##### **Risk of injury**

- ◆ Ensure correct positioning of the patient. This applies in particular when using the different accessory parts.

## 5.2 Positioning the patient on the table

### 5.2.1 Hygiene

- 1 Check the unit for cleanliness.
- 2 If necessary, cover the tabletop with fresh paper.

### 5.2.2 Radiation protection

- ◆ Apply the necessary radiation protection to the patient.

### 5.2.3 Grid and SID

The grid focus (f0) and the source-image distance (SID) must always be identical.

- ◆ Check the correct grid focus (f0) on the grid label before making acquisitions.

### 5.2.4 Table position

Adjust the patient table approximately to the required working height:

- 1 Check the longitudinal and transverse position of the tabletop.
- 2 Adjust the table height.

(→ Page 101 *Preparation and movements of the patient table*)

### 5.2.5 Positioning very heavy patients

#### CAUTION

Too much weight on the table

**Risk of damage or patient injury**

- ◆ When patients weighing more than 200kg are examined, both the patient and the floating tabletop must be centered over the table base.

#### CAUTION

Damage to the tabletop due to overloading, collision or sharp objects

**Risk of injury if the table top splinters or breaks**

- ◆ Do not exceed the maximum weight allowance for the table.
- ◆ Do not allow the patient to sit at the head end of the table.
- ◆ Call Siemens Healthineers Service and have the tabletop checked immediately if there is a possibility that it might have been damaged.

- 1 Center the tabletop.

- 2 Position the patient with a weight between 200 kg and 300 kg only in the center of the tabletop.

### 5.2.6 Tube unit position

- 1 Position the tube unit.
- 2 Observe messages on the display.  
If the tube unit has been rotated about the vertical axis (RVA), the image cutout may not be performed correctly.
- 3 If necessary, correct the image collimation.

### 5.2.7 Manual SID

In cases when the system cannot determine the SID (for example, with oblique projections), the SID must be measured manually.

- 1 Measure the SID using the tape measure.
- 2 If necessary, correct the image collimation.

### 5.2.8 Detector pulled-out

An organ program with detector inside the patient table has been selected and acquisition parameters have been changed. Afterwards the detector has been pulled-out and pushed-in again by mistake. Then the parameters will be reset to their original values.

- 1 Adjust the detector position before selecting an organ program, if possible.
- 2 Check the acquisition parameters set, before releasing acquisitions.

## 5.3 Positioning the patient at the wall stand

### 5.3.1 Hygiene

- 1 Check the unit for cleanliness.
- 2 If necessary, cover the tabletop with fresh paper.

### 5.3.2 Radiation protection

- Apply the necessary radiation protection to the patient.

### 5.3.3 Wireless detector

- Insert the wireless mobile detector in the detector tray.

### 5.3.4 Grid

- ◆ If necessary, remove or insert the grid.

### 5.3.5 Overhead patient hand grip

- ◆ Swing the overhead patient hand grip to the required position or to the side if it is not needed.

(→ Page 168 *Overhead patient hand grip*)

### 5.3.6 Tube unit position

- 1 Press the SmartPositioning button.

(→ Page 206 *Wireless Remote Control Console*)

– or –

Using the handles bring the X-ray tube unit into position, if necessary into the pull-in range of the tracking control.

- 2 Observe messages on the display.

### 5.3.7 Operating Height

- ◆ Adjust the detector unit on the wall stand approximately to the required operating height by pressing the up/down buttons.

## 5.4 IONTOMAT automatic exposure control

The patient table and the wall stand are each equipped with an IONTOMAT measuring chamber for automatic exposure measurement.

All exposures which are taken on the patient table with inserted detector or on the wall stand can be exposed with the IONTOMAT automatic exposure control. The released exposure is then switched off automatically by the IONTOMAT automatic exposure control (1-point technique).

With pulled-out detector, the automatic exposure control cannot be used, exposure must be manually controlled (2-point or 3-point technique).

When selecting a measuring field, make sure the patient is positioned anatomically correctly and according to the measuring field in order to avoid faulty exposures.

For this purpose, three-field templates are available.

- Three-field templates, which can be introduced into the accessory rails of the multileaf collimator for projecting the IONTOMAT dominants onto the object are available for this purpose.

(→ Page 205 *Three-field template, set*)

### 5.4.1 Performing IONTOMAT automatic exposure control

- 1 Select the measuring field according to the object to be imaged.

When collimated narrowly in width, the system switches automatically to the center measuring chamber.

When collimated narrowly in height, no automatic switch-over of the measuring chamber is done. In this case, a smaller amount of imaging radiation will target the lateral measuring chambers.

- 2 If necessary, switch to the center measuring chamber.

No direct radiation may occur close to the measuring field, because the premature shut-off caused by this would lead to underexposure.

- 3 Pay attention to good collimation.

Take care that the correct dominant is selected before the exposure.

For selection of the measuring chamber refer to the Imaging System Operator Manual.

---

The most important part of the object must be located exactly above the measuring field.

Direct irradiation of the area near the measuring field must be avoided, since this would cause premature cutoff resulting in an underexposure.

---

If the measuring field is out side of the collimated field of view (FoV), it will automatically be deactivated by the system.

---

It is not possible to select the left and/or right chamber if one or both of them are not covered by the collimation. In such situation, only the middle chamber can be selected.

An incorrect selection of measuring chambers by the user or by communication network disturbances leads to unnecessary X-ray exposure.

- ◆ Always select the measurement chambers carefully.
  - ◆ Always check the selection of the measurement chambers before starting the measurement.
  - ◆ If several measuring chambers have been selected simultaneously, check the selection carefully.
- 

### 5.5 Additional Cu filter


There are four possibilities of selecting internal additional Cu filters for the automatic collimator:

- 0.0 mm (no) Cu
- 0.1 mm Cu

- 0.2 mm Cu
- 0.3 mm Cu

The selection of additional filters is specified in the organ program. However, you can change this setting at the imaging station or on the TUI.

### 5.5.1 Selecting additional filters

- ◆ Touch this button on the **Basic** tab on the TUI.  
A popup menu appears where you can select the required filter

The additional filter which is selected is displayed at the imaging station.

When changing the additional filters with respect to the organ program selected, this is marked with a star \* at the imaging station, denoting the change.

– or –

Select the filter at the imaging system.

Refer to the Imaging System Operator Manual for further information.

#### CAUTION

Incorrect filter selection

**Risk of increased radiation dose for the patient**

- ◆ Select the filter carefully.

## 5.6 Light localizer and collimation

Ensure that the procedures described in this chapter are performed only by qualified operating personnel.

Collimation close to the object reduces the scattered radiation and thus improves image quality.

The radiation exposure to the patient is also reduced.

To protect the patient from unnecessary radiation, always check the collimation position with the light localizer.

### 5.6.1 Switching on and off of the light localizer

#### Automatic switching on

The light localizer / line light localizer is switched on automatically if

- the collimation size (in the image receptor plane) is changed in width or in height,
- the tabletop brake is released,
- the tube unit is moved in longitudinal direction or in height,
- the detector unit in the patient table is moved,
- the detector unit on the wall stand is moved in height,
- the top alignment function is activated/deactivated.

The service technician can configure when which function is activated individually or together.

#### CAUTION

Looking into the laser light

##### Eye injury

- ◆ Take care that neither you nor anybody else looks directly into the laser beam.
- ◆ Close the sliding cover over the laser light to protect the eyes of the patient and other people.
- ◆ Note that the collimator can heat up when the light localizer is used for a longer period of time.

To protect the eyes of the patient or any other person, the LASER radiation exit of the LASER line light localizer can be closed with the sliding cover.

(→ Page 67 *Multileaf collimator*)

#### Automatic switching off

The light localizer is switched off, independently of how it was switched on, after a configurable time.

The internal time switch can be set by service: 5 s to 118 s

#### Switching on/off manually

You can switch the light localizer on or off at any time independent of automatic switching on and off.

- ◆ Press the button on the collimator.

## 5.7 Collimation

Collimation can be set in the following manners:

- Initial collimation is defined in the organ program.
- Via ACSS  
(→ Page 137 *Collimation / Automatic Collimation Size Sensing*)
- You can set the collimation on the imaging system virtually in the camera image (Virtual Collimation).  
See the Operating Manual of the imaging system.
- You can set the collimation with the light localizer via the knobs on the multileaf collimator.

### 5.7.1 Collimating manually

- 1 Rotate the diaphragm knob to adjust the format height and width.
- 2 Check illumination of the format and collimate as close as possible to the object.

Collimation close to the object reduces the scattered radiation and thus improves image quality. The radiation exposure to the patient is also reduced.

To protect the patient from unnecessary radiation, always check the collimator position with the light localizer.

### 5.7.2 Collimation / Automatic Collimation Size Sensing

MULTIX Impact C supports ACSS (Automatic Collimation Size Sensing).

ACSS limits the maximum opening of the collimator to the borders of the selected detector in the table and Bucky wall stand. ACSS is indicated by the letters ACSS in the collimator display. In case ACSS is not possible, Manual will be displayed. In this case you have to take care that the X-ray field will not go beyond the border of the FD detector.

The collimation format can be set in the OGP. If the selected OGP requires a new position of the tube and/or the detector, the collimation is set when the detector and the tube have reached the final position stored in the OGP.

Conditions for ACSS:

- All system components (collimator, generator, tube stand, imaging system, detector) are ready for use.
- A detector is inserted in the tray of the table or wall stand (ACSS is not supported for free mobile detectors or free cassette operation).
- The central beam is perpendicular to the solid state detector with a tolerance of not more than  $\pm 3$  degrees.

The SID is between 100 cm and 180 cm.

Make sure the X-ray beam does not extend the dimensions of the detector.

The collimator is not rotated.

The respective icon in the collimator status bar is not visible.

Manual collimation is always active if the above conditions are not fulfilled.

The type of collimation is displayed on the TUI.

- **ACSS** is displayed: Automatic Collimation Size Sensing (ACSS) is active.
- **Manual** is displayed: Manual collimation is active.

Make sure the X-ray beam does not extend the dimensions of the detector.

#### CAUTION

Precision of collimator control is not high enough to limit X-ray field inside detector active area when the system is not in **ACSS** mode.

**Patient receives more radiation than necessary.**

- ◆ In **Manual** mode, check collimator size and make sure that the collimator does not open wider than the detector active area.


If **Manual** is active, the user is responsible for the adjustment of the multileaf collimator.

If **Manual** is active, but the image is very flat (low contrast), auto cropping may not work properly. In this case, disable auto cropping at the imaging system, and use the manual shutter after acquisition to crop the image manually.

### 5.7.3 Restoring the last collimation

This button on the TUI serves for memory recall of the last used format. If you press this button, then the collimation set for the last exposure is restored.

Please note that the memory recall can take place only if the preceding exposure was actually released.

## 5.8 Exposures onto the patient table

(→ Page 139 *Vertical projection onto the patient table*)

(→ Page 140 *Vertical projection onto the pulled-out detector*)

(→ Page 141 *Oblique projection onto the patient table*)

### 5.8.1 Vertical projection onto the patient table

#### Positioning the tube unit automatically

- ◆ Press the SmartPositioning button on the Wireless Remote Control Console.  
The tube unit moves to the position provided in the selected organ program.

#### Centering the tube unit manually

Manual movement of the tube unit to the required position:

- 1 Position the tube unit manually as required for the exposure using the movement buttons on the control panel or with the buttons on the handles.
- 2 Bring the tube unit into the centered mid-position.
- 3 Press the centering button on the tube unit to center the detector to the tube.

#### SID

Adjust SID by tube unit or vertical table movement (optional):

- 1 Press button on the tube unit control panel.

– or –

Actuate the table height adjustment via the foot kick switches at the table base or via the footswitch (option for the elevating patient table).

- 2 If required, switch on the tracking control for SID with the SID tracking button.

#### Light localizer

Switch the light localizer on the multileaf collimator on if required:

- ◆ Press the button, if necessary.


#### Positioning the patient

Push the patient with the floating tabletop into the required position:

- ◆ Actuate the foot kick switch or foot switch.

### 5.8.2 Vertical projection onto the pulled-out detector

Example

#### **CAUTION**

Tube is not directed toward the detector during acquisition.

#### **Unwanted radiation exposure**

- ◆ Use the free exposure mode with care.
- ◆ Check that the tube is active and use the light marker for positioning before releasing X-ray.


The automatic exposure control and the grid are not available for free exposures with the detector tray pulled out.

You must therefore expose in 2-point and 3-point technique.

### Preparations

- ✓ An organ program for an exposure on the pulled-out detector unit (table) is selected.
- ◆ Pull the detector tray out up to the stop.  
The acquisition type is switched over automatically.


Example

### Positioning the tube unit

- ◆ Using the handles, bring the X-ray tube unit to the corresponding stop position above the detector and into vertical beam path (0°):  
Loosen the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.

Example for axis marking on the handle

## SID

Set the SID by tube unit or table vertical movement:

- 1 Press the button on the tube unit control panel.

– or –

Actuate the table height adjustment.

- 2 Measure the SID with the tape measure: Distance up to the detector.

– SID 90 cm +

- 3 Adjust the SID value on the TUI.

## Positioning the patient

Push the patient with the floating tabletop into the required position:

- ◆ Actuate the foot kick switch or foot switch.

### 5.8.3 Oblique projection onto the patient table

- 1 Rotate tube unit about the horizontal axis (RHA) and/or about the vertical axis (RVA) into the required direction.
- 2 Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.

The angle of the rotation about the horizontal axis (RHA) is displayed in degrees in the display of the TUI.

Pay special attention in the diagnosis to distortions of the image geometry due to the oblique projection!

Automatic Collimation Size Sensing (ACSS) is not effective in the oblique projection.  
You must collimate manually.

Make sure the collimated field is covered by the detector in the patient table.

3 Make the exposure.

Pay special attention in the diagnosis to distortions of the image geometry due to the oblique projection!

#### CAUTION

Distortion due to exposure geometry

#### **Incorrect basis for diagnosis**

- ◆ Pay special attention in the diagnosis to distortions of the image geometry when using oblique projections.

## 5.9 Exposures onto the wall stand

(→ Page 142 *Horizontal projection onto the wall stand*)

(→ Page 143 *Vertical projection onto the wall stand*)

(→ Page 143 *Oblique projection onto the wall stand*)

### 5.9.1 Horizontal projection onto the wall stand

#### Positioning the tube unit automatically

- ◆ Press the SmartPositioning button on the Wireless Remote Control Console.  
The tube unit moves to the position provided in the selected organ program.

#### Centering the tube unit manually

Manual movement of the tube unit to the required position:

- 1 Position the tube unit manually as required for the exposure using the movement buttons on the control panel or with the buttons on the handles.

Please take care especially for exposures on the wall stand to exact fixation of the overhead support in one of the predetermined SID stop positions.

You avoid image quality defects due to grid shadows if the set SID agrees with the focusing of the grid.

- 2 Center the tube unit to the detector.

(→ Page 115 *Centering the tube unit in vertical direction*)

The tube unit automatically follows the detector in height adjustment within a certain range (tracking control).

## 5.9.2 Vertical projection onto the wall stand

### Positioning the tube unit automatically

- ◆ Press the SmartPositioning button on the Wireless Remote Control Console.  
The tube unit moves to the position provided in the selected organ program.

### Centering the tube unit manually

Manual movement of the tube unit to the required position:

- 1 Position the tube unit manually as required for the exposure using the movement buttons on the control panel or with the buttons on the handles.
- 2 Set the SID by using the buttons on the tube unit control panel.

– or –

Set the SID by using the buttons on the detector unit.

SID tracking is possible for this mode.

(→ Page 99 *Activating SID tracking in wall acquisition mode*)

## 5.9.3 Oblique projection onto the wall stand


### CAUTION

Tube is not directed toward the detector during acquisition.

#### Unwanted radiation exposure

- ◆ Use the free exposure mode with care.
- ◆ Check that the tube is active and use the light marker for positioning before releasing X-ray.

- 1 Tilt the detector unit into the horizontal position.  
The mechanism locks automatically into place.

- 2 Rotate tube unit about the horizontal axis (RHA) and/or about the vertical axis (RVA) into the required direction.
- 3 Press the buttons for the longitudinal, transverse, and vertical brake or for the 'floating' movement and bring the tube unit into the required position.
- 4 Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.

The angle of the rotation about the horizontal axis (RHA) is displayed in degrees in the display of the TUI.

Pay special attention in the diagnosis to distortions of the image geometry due to the oblique projection!

Automatic Collimation Size Sensing (ACSS) is not effective in the oblique projection. You must collimate manually.

- 5 Make the exposure.

Pay special attention in the diagnosis to distortions of the image geometry due to the oblique projection!

#### CAUTION

Distortion due to exposure geometry

##### **Incorrect basis for diagnosis**

- ◆ Pay special attention in the diagnosis to distortions of the image geometry when using oblique projections.

## 5.10 Free exposure onto the wireless detector

#### CAUTION

Tube is not directed toward the detector during acquisition.

##### **Unwanted radiation exposure**

- ◆ Use the free exposure mode with care.
- ◆ Check that the tube is active and use the light marker for positioning before releasing X-ray.

**CAUTION**

Direct radiation on the detector

**Ghost image appears on next image.**

- ◆ Avoid direct radiation on the detector.
- ◆ Wait approximately 30 seconds after such an acquisition before making the next one.

For free exposures onto detector, Automatic Collimation Size Sensing (ACSS) and automatic exposure control are not available.

You must therefore collimate manually and expose in 2-point or 3-point technique. If necessary, the clip-on grid must be attached to the detector.

### 5.10.1 Preparations

- 1 Select an organ program on the imaging system for a free exposure onto wireless detector or for an exposure onto the bed.
- 2 Check which detector is selected on the imaging system. Make sure that the serial number of the wireless detector to be used is identical to the one displayed on the imaging system.
- 3 Make sure that the indication status light on detector is green and so ready for radiation.

### 5.10.2 Positioning the tube unit

Using the handles, bring the X-ray tube unit into position:

- ◆ Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.
  - (→ Page 146 *Positioning for vertical projection*)
  - (→ Page 146 *Positioning for oblique projection*)
  - (→ Page 146 *Horizontal projection onto the mobile detector*)


### 5.10.3 SID

Set the SID by tube unit or table vertical movement:

- 1 Press the button on the tube unit control panel.  
– or –  
Actuate the table height adjustment.
- 2 Measure the SID with the tape measure: Distance up to the detector.

- 3 Adjust the SID value on the TUI.

#### 5.10.4 Positioning for vertical projection

- 1 If necessary set the tube unit vertical.
- 2 Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.

#### 5.10.5 Positioning for oblique projection

Rotate the tube unit about the horizontal axis (RHA) and about the vertical axis (RVA) into the required direction:

- ◆ Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.

The angle of the rotation about the horizontal axis (RHA) is displayed in degrees in the display of the TUI.

Pay special attention in the diagnosis to distortions of the image geometry due to the oblique projection!

We recommend using a protector for wireless detector when patient is directly standing on it.

#### CAUTION

Distortion due to exposure geometry

##### **Incorrect basis for diagnosis**

- ◆ Pay special attention in the diagnosis to distortions of the image geometry when using oblique projections.

#### 5.10.6 Horizontal projection onto the mobile detector

- 1 Attach the detector holder to the required position and insert the detector.  
(→ Page 196 *Detector holder, lateral*)
- 2 Rotate tube unit about the horizontal axis (RHA) and/or about the vertical axis (RVA) into the required direction.
- 3 Release the brakes with the push buttons on the control panel or with the buttons on the handles and position the tube unit.
- 4 Make the exposure.

## 5.11 Pediatric Use: Summary

### 5.11.1 Introduction

Special care should be exercised when imaging patients outside the typical adult size range, especially smaller pediatric patients whose size does not overlap the adult size range (e.g., patients less than 50 kg (110 lb) in weight and 150 cm (59 in) in height, measurements, which approximately correspond to that of an average 12 year old or a 5th percentile U.S. adult female.

Exposure to ionizing radiation is of particular concern in pediatric patients because:

- for certain organs and tumor types, younger patients are more radiosensitive than adults (i.e., the cancer risk per unit dose of ionizing radiation is higher for younger patients);
- use of equipment and exposure settings designed for adults of average size can result in excessive and unnecessary radiation exposure of smaller patients; and
- younger patients have a longer expected lifetime over which the effects of radiation exposure may manifest as cancer.

To help reduce the risk of excessive radiation exposure, you should follow the ALARA (As Low As Reasonably Achievable) principle and seek to reduce radiation dose to only the amount necessary to obtain images that are adequate clinically.

### 5.11.2 Use of anti-scatter grid

The grid can be removed from the side of the detector housing with an easy de-lock mechanism. The information which grid is inserted or if the grid is removed, is shown at the imaging system as well as at the touch user interface at the tube head. It is possible to program if a grid shall be available or not in the organ program, this means the user gets a message to remove the grid for pediatric examinations.

#### CAUTION

Image made with incorrect grid status during pediatric examination

#### Radiation not as desired

- ◆ The expected grid status (inserted, removed) is indicated at the imaging system.
- ◆ Manually remove the grid as required before making the X-ray exposure.

### 5.11.3 Use of filtration

The X-ray beam quality can be changed to reduce low-energy quanta and thus reduce the skin dose. The collimator is equipped with copper filters of 0.1, 0.2 or 0.3 mm Cu. The usage of the filters can be programmed in the organ program. Dedicated pediatric organ programs include this filtering, the filters will be inserted automatically to match the programmed value. The filters can be switched on and off any time at the imaging system as well as at the touch user interface at the tube head.

### 5.11.4 Automatic Exposure Control (AEC)

The system is equipped with a multi-field measuring chamber with sensitivity suitable for pediatric imaging in examinations with automatic exposure control. The different fields can be selected or programmed for the position or anatomy of the current task. The status of the selected chamber is displayed at the Imaging system as well as at the touch user interface.

### 5.11.5 Specific organ programs for pediatric examinations

Our system provides a comprehensive set of dedicated pediatric organ programs. It is possible to add organ programs or change the organ program parameters according to the customers needs. If additional specific organ programs for a specific patient age or body region are needed, please contact the Siemens Healthineers application specialist in your country for support.

For further information, please see the following websites:

- Right Dose Information Center

Find training materials and education across all modalities and clinical specialties.

<https://www.healthcare.siemens.com/medical-imaging/low-dose?stc=wwhim801487>

- Understanding Medical Radiation

A knowledge resource for patients and caregivers.

<http://www.medicalradiation.com/>

Additional educational material such as a checklists that includes special pediatric issues during pre-acquisition, acquisition, and post-acquisition may be obtained from the following example:

"The Alliance for Radiation Safety in Pediatric Imaging, Implementation Manual Image Gently® Digital Radiography Safety Checklist, Safety Steps to Do and Verify for Your Pediatric Patient".

<http://www.imagegently.org/>

### 5.11.6 Pediatric positioning aids

Positioning aids are available that support a safe examination with small patients. These aids are described in Register Accessories and Auxiliary Devices.

## 5.12 Ortho operation

If configured, you can acquire a series of images from different patient heights at the Bucky wall stand or on the patient table.

Both SmartOrtho and Smart Virtual Ortho can provide Ortho acquisition with automatic movement, whereby the system will automatically move in the correct position for each image of the Ortho sequence after the user has selected start and end position.

(→ Page 149 *Ortho on patient table*)

(→ Page 149 *Ortho on wall stand*)

(→ Page 150 Organ program for Ortho)


(→ Page 150 Preparing Ortho)

(→ Page 155 Releasing the Ortho exposures)

### 5.12.1 Ortho on patient table

The Ortho acquisition is done by tilting the X-ray tube unit and moving the detector unit. The X-ray tube unit is centered over the table. When an Ortho sequence is started, the detector moves automatically to the different positions (up to 3). The tube follows the detector, by tilting itself accordingly.

Example: Longitudinal Ortho a.p. lying (region with 3 positions)


Ortho acquisitions on the table can be performed with a SID of 1.5 m.

### 5.12.2 Ortho on wall stand

The Ortho acquisition is done by tilting the X-ray tube unit and moving the detector unit. The X-ray tube unit is positioned in the middle of the Ortho cover range. When an Ortho sequence is started, the detector moves automatically to the different positions (up to 4). The tube follows the detector, by tilting itself accordingly.

Example: Longitudinal Ortho a.p. standing (region with 3 positions)

For automatic Ortho imaging onto the wall stand, the Ortho stand is needed.

(→ Page 171 *Multipurpose / Ortho stand*)

Ortho acquisitions onto the wall stand can be performed with large SID (3 m) and short SID (1.8 m).

### 5.12.3 Organ program for Ortho

Organ programs for Ortho allow up to four acquisitions, belonging to one Ortho series. For each step individual settings for the generator are available.

The position of the units must be saved for the Ortho organ programs. It must be ensured that the unit axes are located orthogonally to each other.

#### CAUTION

Wrong organ program used in Ortho

**Delayed diagnosis or repeated examination because results cannot be composed**

- ◆ Use the correct organ program for Ortho.
- ◆ Make a test run.

### 5.12.4 Preparing Ortho

Pay attention to the messages displayed on the imaging system display and on the TUI.

If the series acquisition is interrupted, the whole sequence must be started again from the beginning!

(→ Page 150 *Organ program for Ortho*)

- ✓ An Ortho organ program is selected.
- ✓ The collimator must be in the 0° lock-in-position.
- ✓ For Ortho acquisitions onto the wall stand, the detector unit must be in the 0° position.

- 1 Press the SmartPositioning button (if configured) to move the X-ray tube assembly and the detector to the start position.

#### CAUTION

The elevating table collides with patient during Ortho positioning.

**Injury to patient.**

- ◆ Ensure that nobody is in the danger zone of the elevating table when using the Wireless Remote Control Console for Ortho positioning.

– or –

If the SmartPositioning function is not configured, position the tube unit in detent and center to the detector.

- 2 Position the patient lying supine on the patient table or standing with their back directed to the Ortho stand.

(→ Page 151 *Positioning the patient on the patient table*)

(→ Page 151 *Positioning the patient on the Ortho stand*)

## Positioning the patient on the patient table

### CAUTION

Patient's hand is under the table during ortho

**Risk of crushing between table and detector tray**

- ◆ Check the position of the patient's hands.
- ◆ Make sure the patient is using the hand grips.

- 1 Position the patient lying on the table.
- 2 Position the calibration element (sphere, ruler) in the area relevant for diagnosis.

## Positioning the patient on the Ortho stand

### WARNING

Use of Ortho functionality without Ortho stand.

**Risk of injury to patient**

- ◆ Always use the Ortho stand if you are doing an Ortho study.

(→ Page 171 *Multipurpose / Ortho stand*)

- 1 Carefully position the Ortho stand as close as possible to the detector tray.

Check that the detector tray does not grind on the semi-transparent stripes of the complete relevant area.

- 2 Position the patient standing on the Ortho stand in such a way that the longitudinal axis of the object runs parallel to the image receptor.
- 3 Position the calibration element (sphere, ruler) in the area relevant for diagnosis.


## Changing/entering the TOD


For later measurements and calculation of acquisition positions, the correct TOD = "Table"-Object Distance must be entered, i.e. Ortho stand to object plane distance.

- or —

- ✓ SmartOrtho is configured.
- ✓ (→ Page 150 *Preparing Ortho*)

| Icon                                                                                | Short description                                                      |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Select TOD (Table-Object Distance)                                     |
|  | Acquisition start or end position (based on the acquisition direction) |
|  | Acquisition end or start position (based on the acquisition direction) |

- 

- 5 Press the other icon on the TUI.

The position is stored and the icon is marked with a check mark.

The acquisition order has been set in the Ortho organ program.

You can change the acquisition order in the **OGP Editor**.

If the upper or lower limit is out of the travel range of the detector, the icon is marked with a cross mark.

- Redefine the ROI.
- Use the additional removable platform on the Ortho stand when performing Ortho acquisition on the wall stand.

(→ Page 171 *Multipurpose / Ortho stand*)

The system automatically calculates any intermediate acquisition positions required for coverage of the acquisition range.

If more than four acquisition positions would be required to cover the acquisition range, a message is displayed. The target area has to be adjusted and the start and end positions have to be set again.

- 6 Turn off the laser light.

- 7 Collimate the region of interest in width.

(Collimation in height will be adjusted automatically.)

- 8 Check acquisition parameters, collimation, and measuring fields.

When collimated narrowly in width, the system switches automatically to the center measuring field.

## Adjusting the ROI with camera in the control room

Using Smart Virtual Ortho with the integrated 2D camera, you can smoothen the preparation for orthopedic procedures, especially long leg and full spine, and thus enhance the conditions that patients experience.

- ✓ Smart Virtual Ortho is configured.
- ✓ (→ Page 150 *Preparing Ortho*)
- ✓ The camera is configured for Ortho.
- ✓ The X-ray tube unit is positioned for Ortho.
- ✓ TOD has been set correctly.

- 1 Click the button on the imaging system.

The Ortho acquisition range (ROI) is shown with an orange frame on the imaging system.

- 2 Select the start and end position by adjusting the upper and lower line.

If the ROI is out of the travel range of the detector, the upper and lower lines are marked in red.

- Redefine the start or end position.
- Use the additional removable platform on the Ortho stand when performing Ortho acquisition on the wall stand.

(→ Page 171 *Multipurpose / Ortho stand*)

- 3 Collimate the region of interest (ROI) in width by adjusting the right and left line.
- 4 Check acquisition parameters and measuring fields on the imaging system.

#### ROI auto detection

- ✓ Smart Virtual Ortho is configured.
- ✓ Auto Long-Leg/Full-Spine Collimation is configured.
- ✓ (→ Page 150 *Preparing Ortho*)
- ✓ The camera is configured for Ortho.
- ✓ Workplace wall stand is selected.
- ✓ The X-ray tube unit is positioned for Ortho.
- ✓ TOD has been set correctly.

- ◆ Click the button on the imaging system.

The ROI is automatically detected and shown with a blue frame on the imaging system.

If the ROI is out of the travel range of the detector, the upper and lower lines are marked in red.

- You can manually adjust the ROI.
- Use the additional removable platform on the Ortho stand when performing Ortho acquisition on the wall stand.

(→ Page 171 *Multipurpose / Ortho stand*)

- Adjust the upper or lower offset in the **Ortho Setting** function.

It is the responsibility of the operator to ensure correct acquisition range. Adjust the ROI manually if necessary.

To adjust the ROI manually:

- Select the start and end position by adjusting the upper and lower line.
- Collimate the region of interest (ROI) in width by adjusting the right and left line.

Once the ROI has been adjusted manually, the frame turns into orange.

### 5.12.5 Releasing the Ortho exposures

- ✓ The first Ortho step is selected.
  - ✓ The units and the patient are correctly positioned.
  - ✓ The Ortho acquisition range (ROI) has been defined.
  - ✓ The TOD and the acquisition parameters are correctly set.
- 1 Confirm all the parameters by clicking **Confirm** on the imaging system.

While using the Ortho function, the X-ray exposure should be performed only after the suggested image field has been validated by the system operator.

To reselect the Ortho organ program, click **End of Exam**.

- 2 Press the acquisition push-button half way (pre-contact) and hold it in that position.  
The system moves to the first acquisition position.

The acquisition order has been set in the Ortho organ program.  
You can change the acquisition order in the **OGP Editor**.

- 3 Instruct the patient to keep still.
- 4 Press the acquisition push-button fully down and keep it pressed until all the acquisitions are finished.  
The first position is acquired and displayed.  
The system moves to the next acquisition position.

The next position is acquired and displayed.

and so on

After all exposures have been carried out, the acquired images are automatically stitched together.

Please pay attention to the system messages displayed on the monitor. These messages indicate the progress of the Ortho acquisitions.

If an Ortho series acquisition is interrupted, the whole sequence must be started again from the beginning!

The same Ortho organ program will be added in the Ortho examination mode automatically.

In this case, the acquisition range and modifications on the organ program parameters will be kept.

However, there is no frame on the imaging system to indicate the acquisition range anymore.

Start acquisitions when the system is ready.

---

## Accepting or rejecting the acquired images

After acquisition, you can accept or reject the image:

- ◆ If the image is okay: Click **Accept**.

– or –

If the image is not okay: Click **Abort**.

The same Ortho organ program will be added automatically.

In this case, the default organ program parameters will be used. Any modifications will be discarded.

---

## 6 Accessories and Auxiliary Devices

### 6.1 Preliminary remarks

#### 6.1.1 Equipment with accessories

Depending on the system type and configuration, the accessories described in the following chapters either from part of the standard equipment supplied or can be purchased optionally.

For more information, please contact your local sales representative.

Please note that not all accessories described in this Operator Manual are deliverable in all countries. For detailed information, please contact your local sales representative.

#### Other products/components

For operation, technical description, models, and technical data, please refer to the documentation supplied by the respective manufacturer.

#### 6.1.2 Proper use of the product

Proper use of the accessories presupposes that the user has read and understood this part of the Operator Manual in detail.

In addition, the user must be thoroughly familiar with the rest of the Operator Manual.

Finally, detailed knowledge of the applicable safety regulations for the medical equipment being used is indispensable. These regulations are listed in (→ Page 14 *System Safety*).

Explanation of "proper use" consists of illustrated descriptions of attaching, using and removing the accessory component.

#### CAUTION

Accessories not firmly attached or not used - Accessories overloaded

#### Patient injury

- ◆ Make sure that accessory parts are attached firmly by pulling and pushing them.
- ◆ Follow the instructions in the Operator Manual.
- ◆ Observe the information regarding storage and maximum weight when using the accessories.

### 6.1.3 Safety

#### WARNING

Wear or material fatigue in the system or accessories

##### **Risk of injury or system damage**

- ◆ Follow the maintenance guidelines to maintain safety and proper function of the system.
- ◆ Check accessories for wear before use.

#### CAUTION

Use of accessories not intended for the system

##### **Risk of injury to patient or user - Risk of damage to the system**

- ◆ Use only accessories approved by Siemens Healthineers with this system.
- ◆ All accessories are labeled and described in the Operator Manual.

For reasons of safety for patients, operators and the equipment, only original accessories manufactured and sold by Siemens Healthineers or Siemens Healthineers-approved accessories from another manufacturer may be used.

Siemens Healthineers cannot assume any responsibility for accessories not approved by us.

#### CAUTION

Trolley or wheeled accessory tips over or suddenly stops while being rolled

##### **Injury of persons**

- ◆ Please be careful when pushing trolleys or wheeled accessories. They are not designed to be pushed over a step, cable or other hindrance because they can tip over or suddenly stop.
- ◆ Take special care when wheeling the Babix cradle trolley with a baby in the Babix cradle.

#### CAUTION

Lead radiation protection can drop on the patient

##### **Risk of patient injury**

- ◆ Avoid any collisions with the radiation protection when moving the system.
- ◆ Make sure that the radiation protection is securely fastened.

### 6.1.4 Orientation

When using the accessories and to simplify the description, the following directional indications shall apply in this Operator Manual with regard to the working position of the user.

- In front means in front of the patient table looking in the direction of the X-ray tube
- Left is the head end of the patient table
- Right is the foot end of the patient table
- Frontal means the longitudinal side of the patient table close to the operating panel
- In back means the longitudinal side of the patient table away from the operating panel
- Turn toward the right means tighten accessory with knurled screws or levers (clockwise)
- Turn toward the left means remove accessory with knurled screws or levers or other fastening devices (counterclockwise)

#### Example

- 1 Turn both knurled screws toward the left.
- 2 Push the grip protection bar on the head end into both accessory rails.
- 3 Turn both knurled screws toward the right.
- 4 Check that the grip protection bar is seated securely.


## 6.2 Description

### 6.2.1 Trolley UM mobile patient table

#### Intended use

As a patient supporting device of the X-ray imaging system, Trolley UM mobile patient table can facilitate X-ray imaging of the recumbent patient, for instance the examination of stomach and lower limbs.

## Product description

- (1) Accessory rails
- (2) Castors without brakes
- (3) Castors with brake pedal

**Castor brake pedal** To enable the patient table to be fixed at any location, the two castors at the front of the table are equipped with brake pedals (as shown in the picture).  
When the patient is moving onto or leaving the patient table, this device can be locked to fix the table position.

- Product features**
- Wide application scope
  - High reliability
  - Easy mobility of the patient table
  - Tabletop has high level of X-ray permeability and stiffness
  - Large tabletop facilitates emplacement of patients
  - The thin tabletop surface can narrow the distance between the cassette and the object to be exposed;
  - By using the brake control device, it is easy to position the patient table.

Trolley UM patient table must not be used to carry patients, it is only to be used as auxiliary supporting device when taking exposures.

When the patient is moving onto or leaving the patient table, the brake pedal must stay locked. Besides, give assistance to the patient, otherwise the table might slide to create hazard.


## Risk of crushing

When the equipment is handled properly and when patients are positioned, the staff may only grasp the handles provided with their hands.

It is the responsibility of the operator to initiate movements of the equipment only if it is ensured that neither the patient nor third parties can be endangered by these movements.

The parts of the unit marked in the drawing below indicate possible danger zones where the patient or staff could crush their fingers or be struck.

Risk of injury for the operating staff when placing the tabletop in position.

## Operation

- Requirements** As a part of the X-ray imaging system, the mobile patient table is used for the examination of body parts such as abdomen and lower limbs.
- The X-ray imaging system must be operated by trained and qualified personnel. In addition, make sure to follow the requirements of relevant laws and regulations during the operating process.

## How to operate the patient table

Before using the table, make sure to remove all nearby objects such as chairs which might collide with the table.

When moving the mobile patient positioning table:

- Make sure to avoid collision with X-ray the detector or other devices.
- Make sure to notice and avoid squeezing of hands and similar incidents.
- Make sure that both the patient and the operating personnel stay inside controllable safety zone.
- Make sure to take care and avoid bringing any unnecessary harm to patients, operators or other personnel, and also tube and detector.

- 1 Before exposure, place the detector and the tube at the required exposure position.
- 2 Adjust the height of the detector moving it at pre-set height matching the patient table height.
- 3 Move the table to a position above the detector, make sure the table cannot collide with the detector and the post.
- 4 Move the table to an appropriate position.
- 5 Lock both castors and ensure the table is fixed.
- 6 Help the patient onto the bed and adjust the position as required.
- 7 Release the brakes of both castors, move the table to obtain the required exposure position, and then lock the brakes again.
- 8 Switch on the light localizer at the collimator, and prepare the final exposure position.
- 9 Adjust the light field of the collimator and release an exposure.
- 10 Release the brakes on both castors and move the table to a position convenient for the patient to leave the table.
- 11 Lock both castors and help the patient from the table.

### WARNING

Grid line artifacts

**Normal grid in used in free exposure mode**

- ◆ Do not use normal grid in the mobile patient table. Use the clip-on grid instead.

## Cleaning and disinfection

Before cleaning and disinfecting the equipment, make sure no liquid can seep into the interior of the equipment.

- Cleaning**
- 1 Do not use water to clean the table, otherwise it might bring rust to the mechanical structure.
  - 2 Do not use corrosive, soluble or abrasive cleaning materials for cleaning.
  - 3 For sprayed paint, plastic and chromium alloy surface, make sure to use clean and dry soft cloth to wipe clean.
- Disinfection**
- 1 Disinfect all components of mobile patient table, including accessories, only with a soft cloth soaked in disinfectant.
  - 2 Do not use spray disinfectant.

### Functional and safety check

- Daily checks**
- 1 Check the unit movement
  - 2 Check the brakes
  - 3 Check the mobility of the castors
  - 4 Check that the necessary devices for patient holding and immobilisation are attached properly to the unit.
  - 5 When the castor brake pedal is in locked status, check if the patient table position is locked.

### Technical description

#### Technical data

|                                 |                                                                                                           |                                              |
|---------------------------------|-----------------------------------------------------------------------------------------------------------|----------------------------------------------|
| Dimensions                      | Tabletop                                                                                                  | 2000 mm ( $\pm 5$ mm) X 660 mm ( $\pm 5$ mm) |
|                                 | Distance from tabletop to the ground                                                                      | 700 mm ( $\pm 10$ mm)                        |
|                                 | Usable width of tabletop between standing posts                                                           | 1900 mm                                      |
| Maximum load of the tabletop    | 135 kg                                                                                                    |                                              |
| Aluminum equivalent of tabletop | < 0.8 mm Al (according to measurement results based on GB 9706.12 (IEC 60601-1-3:2013) (100 kV/3.7 mm Al) |                                              |
| Environmental conditions        | Temperature                                                                                               | +10 °C to +35°C                              |
|                                 | Relative humidity                                                                                         | 30% to 70%                                   |
|                                 | Atmospheric pressure                                                                                      | 700 hPa to 1060 hPa                          |

## 6.2.2 Patient positioning mattress

### Application

The mattress is used for comfortable patient positioning and is attached to the protection strip of the tabletop.


### Radiation protection value

The attenuation equivalent is  $< 0.5 \text{ mm Al}$ .

### Attaching the mattress

- 1 Attach the protection strip at the head end.
- 2 Place the mattress on the tabletop and pull the loops over the hand screws of the protection strip.

#### 6.2.3 Paper roll holder

##### Application

The paper roll holder used to hold paper rolls of different diameters and lengths. It can be attached to the foot end or head end of the tabletop.

##### Attachment

- 1 Push the paper roll holder from the foot end or head end onto the tabletop.
- 2 Tighten/loosen it by turning the screws which can be find underneath the paper roll holder.


#### 6.2.4 Hand grip strip

##### WARNING

Hand grips loosen or are not used

##### Injuries to patient from falling

- ◆ Check that the patient hand grips on the table are in the right position.
- ◆ Ensure that hand grips are tightened prior to use.
- ◆ Make sure that the patient uses the hand grips during system movements or when turning over.
- ◆ Use other positioning aids for immobilizing weak or very unstable patients.


### Application

Hand grip strip which can be attached to front or back accessory rail to assist in positioning the patient.

The maximum permissible load of the hand grip strip is 20 kg.

### Attaching the hand grip strip


#### Example

- ✓ Patient table in horizontal position
- 1 Loosen the right and left cap screws (1) on the hand grip strip.
- 2 Position the hand grip strip from below at the accessory rail as shown in above figure.
- 3 Lift it up and turn it down (see following figure)

We recommend placing the hand grip strip on the far side of the user.


- 4 Tighten both cap screws.
- 5 Make sure that the hand grip strip is secure by pulling and pressing on it.

The hand grip strip provides additional safety for the patient. To prevent injury to the patient's hands when moving the tabletop, make sure that the patient uses the hand grip strip (not the cap screw) and does not hold onto the head or the side edge of the patient table.

**Removal** ✓ Patient table in horizontal position

- 1 Loosen both cap screws on the hand grip strip.

- 2 Turn up the hand grip strip and remove it carefully from the accessory rail.

### 6.2.5 Hand grip, lateral

#### **WARNING**

Hand grips loosen or are not used

#### **Injuries to patient from falling**

- ◆ Check that the patient hand grips on the table are in the right position.
- ◆ Ensure that hand grips are tightened prior to use.
- ◆ Make sure that the patient uses the hand grips during system movements or when turning over.
- ◆ Use other positioning aids for immobilizing weak or very unstable patients.

## Application

The hand grips serve for stable and safe patient positioning.


The hand grips are provided for the safety of the patient.


To avoid injuries to the patient's hands when moving the tabletop, please ensure that the patient uses the grips and does not grasp the ends or sides of the tabletop.

## Attaching the hand grip

The hand grip can be attached to either of the two lateral accessory rails of the patient table.

- 1 Hold the base (1) of the hand grip and loosen it by turning the hand grip (2) counterclockwise.
- 2 Insert the hand grip into the groove of the accessory rail and slide it into the desired position.


- 3 Tighten the hand grip by turning it clockwise.

Please make sure that the hand grips are securely fastened.

Removing the hand grip ♦ Proceed in reverse order.

### 6.2.6 Overhead patient hand grip

#### Application

The overhead patient handgrip serves to stretch the patient during chest examinations.

Insert the overhead patient handgrip on the left or right side of the Bucky unit/ detector tray. You can insert it 0° parallel to the detector tray or 90° vertically in front of the detector tray. The overhead patient handgrip engages in the 0° and 90° positions.

#### **⚠ CAUTION**

Overhead patient handgrip overloaded because patient hanging on it.

#### **Risk of damage to wall stand and injury of the patient**

- ♦ The overhead patient handgrip has no locking mechanism to prevent it from being pulled out from the fixture.
- ♦ Please pay special attention when patients use the grip!
- ♦ The maximum load on the overhead patient handgrip must not exceed the stated value (see label).

**CAUTION**

Patient exerts too much weight on the overhead patient hand grip on the wall stand.

**Risk of injury to the patient and damage to the wall stand**

- ◆ Do not let the patient hang on the overhead patient hand grip.
- ◆ Do not exert more than 25kg weight on the overhead patient hand grips.

**Inserting**

- 1 Move the detector tray until its upper edge is approx. at eye level.
- 2 Push the shaft of the overhead hand grip straight down into the holder on the left or right at the top of the detector tray until the shaft is completely inside the duct.

In this position, you cannot swivel it about its vertical axis.

**Moving and removing****CAUTION**

Overhead hand grip on the wall stand swivels unexpectedly

**Injury of patient**

- ◆ Check that the overhead hand grip is correctly locked and that no red mark is visible.


**CAUTION**

Overhead patient hand grip swivels unexpectedly

**Injury of the patient or user**

- ◆ Check that the overhead hand grip (stretch grip) is firmly attached and that no red mark is visible.
  - ◆ Check that it cannot swivel before allowing the patient to use it.
  - ◆ The maximum load on the overhead patient handgrip is 25 kg.
- ◆ When moving or removing the overhead patient hand grip, pull it vertically upwards out of its holder.

### 6.2.7 BABIX holder

- (1) BABIX holder  
(2) Support arm

The BABIX holder is mounted with the support arm on the wall stand column and can be parked either on the left or on the right.

#### Application

The BABIX holder is used to fasten the BABIX during pediatric examinations of small children. To park the holder, it can be swivelled 180° around its fastening point. The fastening hook slides lengthwise in the support arm. This makes it possible to position the BABIX at different locations up to right in front of the front panel of the wall stand.

#### **⚠ WARNING**

Assembly of babix holder not correctly done.

#### **Injury of persons**

- ◆ Follow the instructions in the Operator Manual.

#### Using the BABIX holder


- 1 Swing the support arm (2) of the BABIX holder (1) in front of the center of the detector tray.

The BABIX holder is now ready to hook up the BABIX.

- 2 Check the BABIX holder (1) for firm attachment at the support arm (2) by pulling or pushing it.

- 3 Check the tractive force as well as screws and other fixing devices of the BABIX holder (1) at the support arm (2) for their stability and loading capacity.
- 4 Hook the bow of the BABIX in the hook of the BABIX holder.

- (1) Attach the BABIX holder on the BWS
- (2) They way for fixing child  
Compression belt can be used to assist fixing child.

#### **WARNING**

BABIX holder loaded with more than 10 kg

#### **Injury of persons**

- ◆ Do not load the BABIX holder with more than 10 kg.

### 6.2.8 Multipurpose / Ortho stand

The Ortho stand is necessary to position the patient correctly for Ortho imaging at the wall stand.

#### **WARNING**

Use of Ortho functionality without Ortho stand.

#### **Risk of injury to patient**

- ◆ Always use the Ortho stand if you are doing an Ortho study.

#### **CAUTION**

Riding on wheeled accessories

**Accessory overbalances and tips over causing damage to system or injury of persons.**

- ◆ Do not ride or allow anyone else to stand, sit or ride on wheeled accessories while they are being rolled.

**⚠ CAUTION**

Positioning the Ortho stand

**Injury of the user**


- ◆ Follow the instructions in the Operator Manual.
- ◆ Be careful not to pinch your fingers, hand or legs during positioning.

**⚠ CAUTION**

Collision of the detector with the Ortho stand


**Risk of damage to the detector and injury of the patient.**

- ◆ Re-position the detector only when no patient is on the Ortho stand.

**Description of the Ortho stand**

- (1) Radiotransparent rear plate
- (2) Accessory rails (on both sides)  
for attachment and height adjustment of accessories
- (3) Hand grips (on both sides) with scales in cm and inch
- (4) Additional removable platform for smaller patients (for example children) and for examinations of the lower extremities
- (5) Brake pedals (on both sides)
- (6) Patient platform (with positioning indentations for additional platform and swivel stool)

## Handgrips on the Ortho stand

- The hand grips (3) are intended to support the patient while standing during the examination.
- The hand grips (3) can also be used as patient stretch grip on top.
- The hand grips can be moved along the guide rails to the required height.
- The hand grips can be used to measure the TOD in cm or inch.
- The maximum load of the hand grips is 30 kg.

### **⚠ WARNING**

Hand grips loosen or are not used

#### **Injury of patient due to falling**

- ◆ Check that the patient hand grips are in the right position.
- ◆ Ensure that hand grips are tightened prior to use.
- ◆ Make sure that the patient uses the hand grips during system movements or when turning over.
- ◆ Use other positioning aids for immobilizing weak or very unstable patients.

### Attaching the handgrips

The hand grips can be attached to the guide rails from below for use as patient support or from above for use as patient stretch grip.

### **⚠ CAUTION**

Mounting the Ortho stand

#### **Patient can slip or fall down**


- ◆ Position the patient hand grips at the lowest possible position when the patient is mounting or dismounting.

**CAUTION**

Patient uses handgrips for support when mounting the Ortho stand.

**Injury of patient**

- ◆ The handgrips can only support 30 kg. If necessary, have the patient grasp the vertical plate when mounting the stand.
- ◆ Then position the handgrips to the correct height for the examination.
- ◆ Make sure that the patient's hands are not on the vertical plate during the examination.


- 1 Turn the locking screws (A) on the handgrips counterclockwise.
- 2 Slide the handgrips from below or from above on the guide rail to the required height.
- 3 Fix the handgrips (A) by turning the locking screws clockwise.

Check that the locking screws on both handgrips are properly fixed. Thus, you prevent the handgrips from disengaging while the patient is leaning on them.

Disengaging handgrips might cause injuries to the patient and/or damage to the support!


In addition to the locking screws, there is a security lever (B) at the bottom of each guide rail to prevent the hand grips from falling to the ground if the locking screws are not properly fixed.

**Removing the handgrips**

- 1 To remove the hand grips when used as stretch grip, just turn the locking screws counterclockwise and remove the handgrips upwards from the guide rails.


- 2 To remove the handgrips when used as support for the patient, turn the locking screws counterclockwise and slide the handgrips downward to the lower end of the guide rails where they are stopped by the security lever.

- 3 Press the security lever to the side as shown in above figure and remove the handgrips downwards from the guide rails.

### Brake pedals

The Ortho stand is equipped with brakes to enable the secure stand of the Ortho stand. For transport of the Ortho stand you have to release the brakes. The brakes are operated with the brake pedals on both sides of the patient platform.

To lock the brakes push the pedal down towards the symbol of the closed padlock (1).

To release the brakes push the pedal down towards the symbol of the open padlock (2).

**CAUTION**

Pedal locking the Ortho stand springs down.

**Risk of foot injury**

- ◆ Make sure that there is nothing under the pedals on the Ortho stand when locking the wheels.
- ◆ Note that there is a pedal on both sides of the stand, and both pedals move when the wheels are locked or unlocked.

**Swivel stool****CAUTION**

The swivel stool detaches from the Ortho stand

**Risk of patient injury**



- ◆ Make sure that the swivel stool is attached correctly to the Ortho stand.
- ◆ Push and pull the stool to make sure that it is secure before positioning the patient.

**CAUTION**

Patient falls from the swivel stool

**Risk of patient injury**


- ◆ Make sure that the patient is seated securely on the swivel stool.
- ◆ Position a handgrip so that the patient can hold on.

- (1) Swiveling seating
- (2) Lever for height adjustment / fixing the swivel movement
- (3) Screw to fix swivel stool in a certain position on the rails
- (4) Rails for longitudinal movement
- (5) Fixing bolt for the attachment of the stool on the patient platform


### Ruler with holding device

The holding device is used to position the ruler on the image.


- (A) Locking screw
- (B) Ruler
- (C) Positioning slots
- (D) Movable holding arm
- (2) Guide rail


#### How to use the holder for the ruler

- 1 Attach the holding device to one of the guide rails (2) and fix it by turning the locking screw (A) clockwise.
- 2 To move the holder vertically, release the locking screw (A).

- 3 The holding arm (D) can be moved horizontally as shown above.


- 4 To position the ruler in one of the two holding slots first of all open the fixing clamp as shown above left.
- 5 Position the ruler in the slot and fix it with the fixing clamp as shown above right.

### BABIX holder

The BABIX holder is used to fasten the BABIX during pediatric examinations of small children.

- (1) Cap screw
- (2) Support arm
- (3) BABIX holder
- (4) Guide rail

#### Attachement

#### WARNING

Assembly of babix holder not correctly done.


#### Injury of persons

- ◆ Follow the instructions in the Operator Manual.

The BABIX holder can only be mounted on the right guide rail at the top.

- ◆ Slide the support arm (2) with the BABIX holder (3) from above on the guide rail (4) and fix it with the cap screw (1).

## Positioning the BABIX holder

- (1) Guide rail
- (2) BABIX holder
- (3) Hook for attaching the BABIX
- (4) Holding knob

## Horizontal positioning

The BABIX holder can be positioned moving it back and forth in the guide rail (2).

## Turning the hook

- ◆ The hook of the BABIX holder can be turned in steps of 90° around the vertical axis by lifting the BABIX holder with the holding knob and turning it.

**⚠ CAUTION**


Babix holder overload

**Risk of damage to of the babix holder and injury of the patient**

- ◆ Do not load the BABIX holder at the point of suspension with more than 10kg.

## Detector holder

- (1) Fixing device for detector
- (2) Detector holder
- (3) Handle
- (4) Locking screw

### Attaching the detector holder

#### CAUTION

The detector holder is not securely fastened.

**Injury of persons or damage to the detector holder**

- ◆ Make sure the detector holder is securely fastened.

The detector holder can be attached only to the right guide rail from below (recommended) or from above.

- 1 Turn the locking screw (4) counterclockwise.
- 2 Slide the detector holder from below or from above on the guide rail to the required height using the handle (3).
- 3 Fix the detector holder by turning the locking screw clockwise.

Check that the locking screw is properly fixed.

In addition to the locking screws, there is a security lever at the bottom of each guide rail to prevent the detector holder from falling to the ground if the locking screw is not properly fixed. (For more details to the security lever, refer to (→ Page 173 *Attaching the handgrips*).

### Placing the detector

- 1 Lift the fixing device (1).
- 2 Place the detector on the lower holder.
- 3 Slide the fixing device downward.
- 4 Ensure the detector is seated properly in the holder.

Removing the detector holder ♦ Refer to (→ Page 173 *Attaching the handgrips*).

### Cleaning of the Ortho stand

#### CAUTION

Patient hand grips slip because the guide rails are contaminated.

##### **Patient injury**

- ♦ Clean the Ortho stand with a suitable agent.
- ♦ Check that there is no oily residue on the guide rails.


### Technical data of the Ortho stand

|                                                                                     |                                                      |
|-------------------------------------------------------------------------------------|------------------------------------------------------|
| Height                                                                              | 197 cm                                               |
| Depth                                                                               | 86 cm                                                |
| Width                                                                               | 82 cm (incl. brake pedals)                           |
| Weight                                                                              | 112 kg (incl. additional platform and two handgrips) |
| Max. patient weight                                                                 | 200 kg                                               |
| Max. patient size                                                                   | 190 cm                                               |
| Attenuation equivalent of the transparent protection plate acc. to DIN EN 60601-1-3 | < 1 mm Al - 100 kV 3.6 HVL                           |

### Location of labels

All labels displayed in the operator manuals are examples only and may differ from the labels attached to the system and components.


(1) Weight label (200 kg)

(2) Type B equipment label

(3) Refer to manual label

(4) Ortho stand ID label

(5) Weight label (120 kg)

This product complies with DHS regulations  
21 CFR Subchapter J, applicable at date of  
manufacture.


Manufactured:

Siemens Shanghai Medical Equipment Ltd.,  
278 Zhou Zhu Road, 201316 Shanghai  
China

(6) DHS label

<1.0 mm Al - 100kV 3.6mm Al HVL

(7) Al equivalent label

(8) Handgrip ID label

(9) Weight label (30 kg)

(10) Caution label


## 6.2.9 Compression Belt

### Safety

#### Risk of crushing

It is the responsibility of the operator to attach and remove the equipment only if it is ensured that neither the patient nor third parties can be endangered.

The parts of the unit marked in the drawing show possible danger zones where fingers can be crushed.

Risk of injury for the personnel when attaching compression belt.

### Functional and safety check

- Daily checks**
- 1 Check tension lever.
  - 2 Check cap screw.
  - 3 When the cap screw and tension levers are in locked status, check if the compression belt is locked.


Remove all traces of contrast agent from the compression belt. Otherwise, the X-ray image shows artifacts.

### Cleaning and disinfection


- Cleaning**
- 1 The compression belt is made of transparent plastic and can be cleaned with a damp cloth.
  - 2 Do not use corrosive, soluble or abrasive cleaning materials for cleaning.
  - 3 For sprayed paint, plastic and chromium alloy surface, make sure to use clean and dry soft cloth to wipe clean.
- Disinfection**
- 1 Disinfect all components of the compression belt only with a soft cloth soaked in disinfectant.
  - 2 Do not use spray disinfectant.

## Product description


- (1) Ratchet handle
- (2) Tension lever
- (3) Cap screw
- (4) Tension lever
- (5) Tag for ratchet handle
- (6) Tag for cap screw

- (1) upward: tightened
- (2) downward: released

Tag for ratchet handle

- (1) upward: released
- (2) downward: tightened

Tag for cap screw

**Intended use** The compression belt is a mechanical component for X-ray exposures.

It is used to quickly fix the patient on the table.

In addition, it can reduce the thickness of the patient which leads to better image quality for some examinations, especially abdominal.

### Product features

- High reliability
- Easy attachment


- Fasten frail patient
- Increase image quality

## Operation


**Requirements** The X-ray imaging system must be operated by trained and qualified personnel. In addition, make sure to follow the requirements of relevant laws and regulations during the operating process.

### For use on the table

- Attachment**
- 1 Move the patient table to the horizontal position.
  - 2 Position the patient on the table.

- 3 Push the tensioning part of the compression belt into the required position over the front accessory rail.


- 4 Press both tension levers (2) outward.
- 5 Check that the tensioning roller is firmly seated by pulling on it.
- 6 Press the ratchet handle (1) downward.

The belt is released.


For obese patients, the belt holder is easier to secure if you are standing at the back of the patient table. In this case, and for restless or frail patients, we recommend that someone assist you in attaching the compression belt

- 7 Lift the holder with the belt across the patient.

- 8** Push the holder over the back accessory rail. Make sure that the compression belt attaches straight across the tabletop and not at an angle.


- 9** Press the tension levers (4) outward.  
**10** Check that the holder is seated firmly by pulling on it.  
**11** Turn the cap screws (3) on the tensioning roller toward the left to tighten the belt.  
**12** Push and pull back on the ratchet handle (1), until the belt tension is optimal.  
**13** Check the belt tension on the patient once again.

#### Loosening the compression belt


- 1** Move the patient table to the horizontal position.
- 2** Press the ratchet handle (1) back.
- 3** Turn the cap screw (3) on the tensioning roller toward the right.  
The compression belt unclenches.
- 4** Use the ratchet handle to set a lower belt tension.

#### Removal

- 1** Move the patient table to the horizontal position.
- 2** Press the ratchet handle (1) back.
- 3** Turn the cap screw (3) on the tensioning roller toward the right.  
The compression belt unclenches.

- 4 Press both tension levers (4) on the belt holder inward.
- 5 Pull the belt holder backward off the accessory rail.
- 6 Lift the belt holder over the patient.
- 7 Roll up the compression belt halfway using the cap screw.
- 8 Place the belt holder next to the tensioning roller.


- 9 Press both levers (2) inward on the tensioning roller.
- 10 Pull the tensioning roller from the front accessory rail.
- 11 Help the patient to get off the patient table.
- 12 Use the cap screw (3) to roll up the belt on the tensioning roller.

### Technical data

| Dimensions |          |
|------------|----------|
| Width      | ≤ 100 mm |
| Length     | ≤ 400 mm |
| Height     | ≤ 150 mm |

| Ambient conditions (Operation)             |                                                   |
|--------------------------------------------|---------------------------------------------------|
| Temperature:                               | +10°C to +35°C                                    |
| Relative humidity                          | 20% to 75% rel. air humidity without condensation |
| Barometric pressure                        | 700 hPa to 1060 hPa                               |
| Ambient conditions (Transport and storage) |                                                   |
| Temperature                                | -20°C to +70°C                                    |
| Relative humidity                          | 10% to 95% rel. air humidity without condensation |
| Barometric pressure                        | 500 hPa to 1060 hPa                               |
| Aluminum equivalent                        | ≤ 0.8 mm Al                                       |

### Location of labels

(1) Identification label for compression belt

(2) Warning label

(3) UR label/ UR


(4) Type B equipment label

### Complying standards

- UL 60601-1, 1st Edition, 2006-04-26 (Medical Electrical Equipment, Part 1: General Requirements for Safety)
- CAN/CSA-C22.2 No. 601.1-M90, 2005 (Medical Electrical Equipment, Part 1: General Requirements for Safety)
- CAN/CSA-C22.2 No. 60601-1, Medical Electrical Equipment - Part 1: General Requirements for Basic Safety and Essential Performance
- ANSI/AAMI ES60601-1 (2005 + C1:09 + A2:10), Medical Electrical Equipment, Part 1: General Requirements for Basic Safety and Essential Performance
- IEC 60601-1:1988 + A1:1991 + A2:1995
- IEC 60601-1: 2005 + A1:2012

### 6.2.10 Cassette holder, lateral

#### Application

- For lateral exposures with a second X-ray tube;
- To hold the cassette upright in all table positions, including Trendelenburg.

#### Maximum cassette format

35 cm in height, width is unrestricted, for example 35 cm x 43 cm horizontal is possible.

#### Minimum cassette formats


13 x 18 cm vertical or 18 cm x 24 cm horizontal.

#### Maximum load capacity

3 kg

#### Attachment

- ✓ Patient table in horizontal position.
- 1 Position the patient.

- 2 Open the clamping lever (1).
- 3 Place the cassette holder on the accessory rail and move it to the exposure position.
- 4 Close the clamping lever (1).
- 5 Push and pull the cassette holder to ensure that it is seated properly.
- 6 Lift the tensioning device (4).
- 7 Place the cassette on the lower holder.
- 8 Slide the tensioning device downward.
- 9 Ensure the cassette is seated properly in the holder.
- 10 Loosen the cap screw (2) to permit adjustment of the height and the cap screw (3) to permit lateral adjustment.
- 11 With the patient positioned for the examination, adjust the height and lateral tilt of the cassette.
- 12 Tighten the cap screws (2) and (3).
- 13 Ensure that the patient is stabilized and comfortable.


Have the patient place the hands on the protection strip, the handgrip strip, or the handgrip, depending on the type of examination. This will increase the feeling of security of the patient during the examination.

## Removal

- ✓ Patient table in horizontal position.
- 1 Hold the rail for height adjustment.
- 2 Loosen the cap screw (2).
- 3 Lift the rail for height adjustment.
- 4 Tighten the cap screw (2).
- 5 Open the clamping lever (1).
- 6 Remove the cassette holder from the accessory rail.

- 7 Assist the patient in stepping down from the patient table.

### 6.2.11 Cassette holder without retainer


#### Application

- For lateral exposures;
- To hold the cassette upright;
- For all cassette formats, horizontal and vertical

#### Attaching the cassette holder without retainer

- 1 Position the patient.
- 2 Ensure the patient is comfortable and stabilized.

- 3 Slide the cassette into the bracket of the cassette holder.
- 4 Ensure that the cassette is seated properly in the holder.
- 5 Position the cassette holder on the tabletop at the patient in the location appropriate for the examination.
- 6 Ensure that the cassette holder is stable.

Have the patient place the hands on the protection strip, the hand grip strip, or the hand grip, depending on the type of examination. This will increase the feeling of security of the patient during the examination.

- Removal**
- 1 Remove the cassette from the bracket of the cassette holder.
  - 2 Remove the cassette holder from the tabletop.
  - 3 Assist the patient in stepping down from the tabletop.

### 6.2.12 Mobile detector holder CR

Refer to the original Operator Manual of the manufacturer.

#### CAUTION

Lifting the table when the mobile detector holder is positioned above the tabletop.

#### Risk of crushing

- ◆ Be careful lifting the tabletop when the mobile detector holder is positioned near the table.

#### CAUTION

Removing the detector from the mobile detector holder - Due to the spring loading, the holder will move upwards when the detector is removed or the brake is released.

#### Risk of injury

- ◆ Be careful when removing the detector from the detector holder.
- ◆ Turn the knob to lock the vertical position before removing the detector.
- ◆ Hold the horizontal arm while slowly removing the detector to make sure that the holder remains at that position.

#### CAUTION


The detector holder is not securely fastened.

#### Injury of persons or damage to the detector or holder

- ◆ Make sure the detector holder is securely fastened.


### 6.2.13 Detector holder, lateral

#### Application


Lateral exposures with the X-ray tube onto the upright detector.

Examples of application

The detector holder is to be used only on or above the tabletop or directly at the side of the tabletop (see examples above).

#### Misuse

- Do not use the detector holder as shown above.
- Do not lean on it, you might damage the tabletop.

#### Attaching detector holder

The detector holder can be attached to either of the two lateral accessory rails of the patient table.

- 1 Take the detector holder with both hands (high weight!) and insert it into the groove of the accessory rail.


- 2 Slide it into the desired position.
- 3 Tighten the detector holder.

**CAUTION**

The detector holder is not securely fastened.

**Injury of persons or damage to the detector or holder**

- ◆ Make sure the detector holder is securely fastened.


- 4 To move the detector holder vertically use the crank (1).

#### Inserting the detector

There are two grooves (A) that can be used to insert the detector:

- In the front groove only the detector can be inserted.
- In the back groove the detector with the optional grid attached can be inserted.


- ◆ Insert the detector into the detector holder from above.

#### Removing detector and detector holder

- 1 Take the detector out of the detector holder from above.
- 2 Place the detector in a safe location.
- 3 Dismantle the detector holder with both hands from the lateral accessory rail of the patient table.

### 6.2.14 Clip-on grids

The following clip-on grids are available:

- 5/85 F115
- 15/80 F115 (only for MAX wi-D)

#### Application

The clip-on grids are used for free exposures on the mobile detectors.



For exposures on the patient table or on the wallstand, use the corresponding internal grid.

When the clip-on grid is not in use, remove it from the detector and store it in a safe place where it cannot fall, for example in the grid holder.

Even if it is not extremely damaged, maybe its characteristics are changed, which may cause a problem in image quality.

#### Attaching the clip-on grids


- (1) Core XL
- (2) MAX wi-D
- (3) Clip-on grid

✓ Detector front must face the backside of the clip-on grid.


- 1 Put the detector with the bottom in the clip-on grid as show in above figure.

- 2 Push the detector gently down against the loading spring of the clip-on grid and let it slide smoothly into the clip-on grid and let go.


Ensure that the grid unit is locked before lifting the detector unit. Otherwise, the grid unit can fall and be damaged.

### Removing the clip-on grid

- ◆ Push the detector gently down against the loading spring of the clip-on grid and then pull the detector out of the clip-on grid and let go.

### 6.2.15 FD protector (patient support)

### Application

The FD protector is a protective cover for the detector. It is used for examinations where the patient has to stand on the detector if the patient body weight exceeds 100 kg.

The FD protector may only be used without clip-on grids.

The FD protector is made of acrylic glass.

The maximum load on the FD Protector must not exceed the specified value on the label.

### Al equivalent

The Al equivalent value is 1.39 (for 3.7 mm Al).

#### 6.2.16 Wall holder for grids

This holder is used for storage of exchangeable grids or cassette trays. You can mount it to the wall.



#### 6.2.17 Compensation filters

##### Application



Absorption compensation for acquisitions of:

- Pelvis
- Foot
- Infant skull
- Adult skull
- T-spine, L-spine
- Shoulder




Compensation filter for pelvis and compensation filter for foot

Compensation filter for infant skull and compensation filter for adult skull

Compensation filter for thoracic spine, lat. and compensation filter for thoracic or lumbar spine, lat.

Compensation filter for shoulder

### Attaching compensation filters


The operating personnel has to take care that the filters are appropriately attached.

- 1 Push the filter with the aluminum part downwards into the two accessory rails of the multileaf collimator.

The safety spring in the left rail is pushed beside in this case.

- 2 Push the filter up to the stop.

**WARNING**

Falling objects because too much weight is hung on the collimator

The rails could be damaged by the use of excessive force. A proper seating of the accessories or auxiliary devices can no longer be guaranteed and can lead to patient injury. The position of the tube assembly can change.


- ◆ Do not hang more than 7 kg on the collimator accessory rails.
- ◆ Do not exert any force on the multileaf collimator or its rails when inserting the accessory or auxiliary device.
- ◆ Check that the accessory is firmly anchored in the collimator rails.

If you exert too strong a pressure on the accessory rails, you might have to change the settings of the RHA because of a change in the position of the tube assembly.

### Removal 1 Pull filter out of the accessory rails of the multileaf collimator

The safety spring in the left rail has to be pushed beside in this case.


### 2 Keep filter secured in the holder.

Please be careful when handling the compensation filters. They are thin, scratch-sensitive, and may become unusable due to rough handling.

### 6.2.18 Filter holder for eight filters

## Application

**For storing a maximum of eight compensation filters.**

### Filter holder operation

### Fastening the wall holder

- ◆ Have the holder fastened to the wall at working height at a suitable place with the enclosed wall plugs and screws.


### Equipping the wall holder

- 1 Remove compensation filters from the packaging.**
- 2 Turn the filter so that you can read the designation on the filters reading normally from above.**
- 3 Place the filters in the upwards angled compartments.**

### Storing the compensation filter

- ◆ Place the filter in a free compartment of the wall holder, turn the filter so that you can read the designation on the filter reading normally from above.

### 6.2.19 Three-field template, set

#### Application

For exposing the IONTOMAT ionization chambers on the object to be exposed.

The three-field templates are available as a complete set or individually for the following SIDs:

- 90 cm - 110 cm
- 110 cm - 130 cm
- 130 cm - 175 cm
- 175 cm - 220 cm

#### Using the three-field template set

- 1 Push the locking lever on the left accessory rail to the left.
- 2 Slide the template into the collimator accessory rail in the correct direction for exposing the ionization chamber.  
The locking lever on the accessory rail springs to the right.
- 3 Check that the template sits firmly in the collimator.
- 4 Expose the ionization chambers. To do so, switch on the light of the collimator.

#### Storing the three-field template set

- 1 Push the locking lever on the left accessory rail to the left.
- 2 Pull the template out of the accessory rails.  
The locking lever on the accessory rail springs to the right.

- 3 Store the three-field templates in a suitable location.

Please be very careful with the three-field templates. They are thin, sensitive to scratching and can become useless if handled carelessly.

## 6.2.20 Wireless Remote Control Console

### Important information

The functions of the Wireless Remote Control Console are activated via an RF transmitter.

The Wireless Remote Control Console may be interfered by other equipment, including portable and fixed RF communication equipment, even if such equipment meets the applicable emission requirements.

The operator must assure that other wireless devices in the 2.4-GHz ISM band are not operated in the vicinity of appr. 5 m around the MULTIX Impact C system. System movements can interrupt sporadically.

Please observe and verify normal operation of the Wireless Remote Control Console prior to using it.

During operation, the Wireless Remote Control Console emits electromagnetic energy. Maybe, nearby electronic equipment is affected.

- Operate the Wireless Remote Control Console only inside the examination room in close vicinity to the tube stand and with full visibility to the patient on the table or at the wall stand.
- Take care that the batteries are always sufficiently charged.

#### WARNING

The Wireless Remote Control Console generates an electromagnetic field so high that it disturbs a life-supporting device.

#### Interference with a life-supporting device resulting in possible malfunction

- ◆ Make sure to keep a safety separation distance of larger than 90 cm between the wireless remote control and the life-supporting device.
- ◆ In case of interference with other equipment increase the distance between the interfering devices.

**⚠ CAUTION**

Collision with operator or patient

**Mechanical or personal damage**

- ◆ Take care that the used remote control and the system have the same colored label.
- ◆ Use remote control only when there is visual and acoustic contact with the patient.
- ◆ In case of unintended movement press the nearest emergency STOP button.

**⚠ WARNING**

Movement without intentional activation of control elements

**Risk of collision, risk of injury to patient or operator - Risk of damage to the system**

- ◆ Do not load the remote console with any objects, accessories, folders or documents.
- ◆ If movements does not stop, press the nearest emergency STOP button.


Changes or modifications not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.

The Wireless Remote Control Console may be interfered with by other equipment, including portable and mobile RF communication equipment, even if this equipment meets the applicable emissions requirements.

In order to perform its intended function, the Wireless Remote Control Console must emit electromagnetic energy. Maybe, nearby electronic equipment is affected.

**Operation**

The remote control is active in a range of approx. 7 m around the tube stand, where the receiver is mounted. If the remote control is in a wider range of the table, it may lose contact with the receiver. It then turns into an inactive mode.

#### Wireless Remote Control Console

- (1) Display area
- (2) Operation area
- (3) Warning label

Example of the warning label

Please ensure that the distance to the next emergency stop button is less than 2 m when pressing a button on the remote control.

**Display area** There are 3 LEDs showing the status of the remote control console.

- LED off: battery sufficiently charged, remote control not connected to the USB charging cable
- LED on red : battery needs charging, remaining capacity < 2.0V
- LED flashing yellow: remote control connected to the USB charging cable and batteries are being charged
- LED on yellow: remote control connected to the USB charging cable and batteries are fully charged


- LED off: The bluetooth signal strength is good.
- LED on: no bluetooth signal

- LED off: No operation button is being pressed or the dead-man sensor is not activated.
- LED on: An operation button is being pressed and meanwhile the dead-man sensor is activated.

## Operation area

(1) Dead-man sensor

On both sides of the remote control, there is a dead-man sensor. To carry out an operation using a button in the operation area, you must activate the dead-man sensor with your hand at the same time.

As soon as you remove your hand from the dead-man sensor, the current operation is stopped.

## Collimator adjustment / tube unit movement

You can select which part of the system you want to operate.

The green LED on the selected button is on.

Selection of collimator adjustment

Selection of tube unit movement

Collimator adjustment is selected by default setting.

Mode-switch can be set by the Siemens Healthineers service to shift the selection from tube unit movement to collimator adjustment after a timeout.

- Opening the collimator in vertical direction

- or -

- Moving the tube unit up in vertical direction

- Closing the collimator in vertical direction

- or -

- Moving the tube unit down in vertical direction

- Closing the collimator in horizontal direction

- Opening the collimator in horizontal direction

### Position of organ program or parking position

SmartPositioning (optional):

Moving the system to the position which is stored in an organ program  
or

Parking position:

Moving the system to the parking position

(The parking position can be configured by Siemens Healthineers service.)

### Detector unit up/down movement

Moving the detector unit of the wall stand up and down

### Tracking control

Disable/enable tube tracking in wall mode

- LED off: Tracking function is inactive.
- LED on: Tracking function is active.

To enable the tracking function, the preconditions must be fully fulfilled.

(→ Page 97 *Tube tracking*)

## Charging

We recommend charging the batteries every day after deployment using the USB charging cable.

USB charging cable

## Changing batteries

The batteries of the Wireless Remote Control Console should be changed every 24 months by authorized and trained staff. Please call Siemens Healthineers Service. Unauthorized exchange may cause severe damage.

## Technical data

|                                    |                                                                                                                                                                                                                                               |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Active range                       | Approximately 7 m around the receiver mounted to the stand.                                                                                                                                                                                   |
| Charging time                      | Approx. 8 h                                                                                                                                                                                                                                   |
| Operating time with full batteries | Approx. 20 h                                                                                                                                                                                                                                  |
| Frequency range                    | 2402 MHz - 2480 MHz                                                                                                                                                                                                                           |
| Modulation                         | GFSK                                                                                                                                                                                                                                          |
| RF Power                           | 4 dBm                                                                                                                                                                                                                                         |
| RF emissions                       | CISPR 11 Compliance: Group 2 FCC 47 CFR 15, RED, SRRC<br>The deployment of Wireless Remote Control in combination with MULTIX Impact C wi-D changes the classification of MULTIX Impact C wi-D from compliance group 1 to compliance group 2. |


## Marking signs

The Wireless Remote Control is marked with a colored sign of different shape. So you can avoid a mix up for installations with a multitude of Wireless Remote Controls. The same mark is attached to the system during installation.

These signs are examples of colored signs of different shape to mark the remote control to the system.

### 6.2.21 Foot switch

The foot switch can only be used with the elevating patient table.

- Left foot switch for table down movement.
- Middle foot switch for table up movement.
- Right foot switch to release the brakes for tabletop floating.


### ⚠ CAUTION

Operator or patient trips over the foot switch cable or the foot switch itself.

#### Risk of injury

- ◆ Position the foot switch and the associated cable on the floor so that nobody trips.


## 6.2.22 Battery charger (MAX wi-D)

- (1) Charger
- (2) Batteries
- (3) 4 charging LEDs for each battery
- (4) Power input with green power indicator  
Green LED is illuminated when power is on.

### Inserting battery into charger

The battery charger is designed to insert a battery in only one direction.

- 1 Insert the battery into the charger as shown above.
- 2 Pay attention to the alignment identification (1) on the battery and the charging slot (see above).

- 3 Gently push the battery onto the connector.

Once the battery is in place correctly, the bottom green LED starts blinking.


Charging a fully discharged battery will take a maximum of 3 hours. Up to three batteries can be charged at a time in the charger.

## Status LEDs

For each charging slot, there are 4 LEDs. When there is no battery to be charged, all LEDs should be off.

There is 1 bottom LED (orange or green) and 3 green LEDs.

Charging levels

Fault signal

## Charging trouble shooting

**ORANGE Fault LED on** The battery is not inserted correctly / battery temperature too high.

- ◆ Remove battery and insert it again correctly.

If the problem remains, the battery may be damaged.

– or –

Remove battery and insert it again with a normal temperature.

The charger slot will remain in a fault condition until the battery is removed and inserted again with a normal temperature.

**No green power LED active** ◆ Check that the power supply is working.


**No charging LEDs active** No battery inserted

Battery not inserted correctly.

- 1 Remove battery and insert it again correctly.

- 2 If problem remains, check that connectors are not damaged or battery is not faulty.


### 6.2.23 Battery charger (Core XL)

- (1) Charger
- (2) Battery
- (3) Power input
- (4) Power input LED  
Green LED is illuminated when power is on.
- (5) Charging LED  
Yellow LED is illuminated when the battery is in charge.
- (6) Battery full LED  
Green LED is illuminated when the battery is full.

#### Inserting battery into charger

The battery charger is designed so that a battery can be inserted in only one direction.

- 1 Insert the battery into the charger as shown above.
- 2 Pay attention to the alignment identification (1) on the battery and the charging slot (Output1) (2) .
- 3 Slide battery into the recess area of the device.  
Once the battery is in place correctly, the bottom green LED starts blinking.

Charging a fully discharged battery will take a maximum of 4 hours.

### 6.2.24 Wall charger

The wall charger has two functions:


- It can store one MAX wi-D portable detector, one MAX mini, one portable grid for the MAX wi-D and a spare detector battery.
- It is also used to charge the battery of the stored MAX wi-D automatically, when the detector is properly inserted in the wall charger. This is indicated by a blinking battery LED on the detector.

The MAX mini and the spare battery cannot be charged.

Do not place the charger in the patient area!

The wall charger can be mounted to the wall.

- (1) Compartment for one MAX wi-D
- (2) Compartment for one MAX mini (without charging)
- (3) Compartment for one spare detector battery (without charging)
- (4) Compartment for one portable grid for the MAX wi-D

#### Inserting the MAX wi-D detector correctly for charging:

- 1 Insert the MAX wi-D detector into the corresponding compartment (charging station) of the wall charger.
- 2 Make sure that the handle of the detector points up and the charging contacts of the detector face the wall to permit the charging procedure to start.

### Removing the MAX wi-D detector:

- ◆ Grab the MAX wi-D detector by its handle and pull it out of the wall charger carefully.

#### 6.2.25 Barcode scanner

The barcode scanner is used for scanning patient data.

The barcode scanner beam must cross the entire label. ( The scanner cannot correctly read a barcode if the entire label is not scanned.)

When connecting the scanner to the imaging system or disconnecting it from the imaging system, make sure the system has been switched off.

### Usage


- ✓ **Pre-registered Patients** list is displayed on the screen.
- 1 Point the barcode scanner at the corresponding label and press the trigger button.  
The patient ID is displayed in the searching field in the **Pre-registered Patients** list.
- 2 Click **Search** .  
The patient is shown in the **Result List**.

**Troubleshooting** If the scanner is not operating properly, the following checks should be performed:

- 1 Check that the interface cable is securely attached to the scanner handle and the terminal.
- 2 Check that the labels are of sufficient quality to be recognized by the scanner. Wrinkled, smudged or torn labels can cause the scanner to not read at all. If you suspect that label quality is a problem, scan a known good label to check the scanner's read operation.
- 3 If the scanner still does not function properly, contact your Siemens Healthineers Customer Service.

## 6.2.26 Intercom system in control room

### Overview

- (1) Loudspeaker
- (2) Microphone
- (3) Listening (loudspeaker) volume (in the control room)
- (4) Speaking (microphone) volume (in the control room)
- (5) LED: Power status indication  
LED lights up when power is ON
- (6) LED: Listening (loudspeaker) status indication  
LED lights up when the loudspeaker is ON
- (7) LED: Speaker (microphone) status indication  
LED lights up when the microphone is ON
- (8) Power ON or OFF button
- (9) Listening (loudspeaker) ON or OFF button
- (10) Speaking (microphone) ON or OFF button



## 7 Technical Description

### 7.1 Location of labels

All labels displayed in the operator manuals are examples only and may differ from the labels attached to the system and components.

#### 7.1.1 Labels on the ceiling-mounted tube support

##### (1) System ID label

The following symbols are contained in the labels, if applicable:


Medical device

## CONFIG

Configurable device. A configurable device is a device that consists of several components which can be assembled by the manufacturer in multiple configurations. Those individual components may be devices in themselves.

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

(6) Power label

(7) China RoHS label

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
(8) Warning label for ionizing radiation for deliveries to China

(9) SGS label

(10) Label for deliveries to Singapore

(11) Brazil label

(12) MoH label for deliveries to Ministry of health in Saudi Arab

(13) EAC label (This label is for Eurasian Economic Union (EEU) with EAC declaration done.)


This product complies with DHHS regulations  
21 CFR Subchapter J, applicable at date of  
manufacture.  
Manufactured August 2018  
Siemens Shanghai Medical Equipment Ltd  
278 Zhen Zhu Road, 201318 Shanghai  
China

(14) DHHS label for deliveries to USA

**MULTIX Impact C**

(15) Malaysia label

### 7.1.2 Labels on the tube unit with multileaf collimator

(1) Refer to manual label


Only for deliveries to China


or

(2) Laser radiation label for collimator AL04 II-D

(3) DAP label for deliveries to Germany and Austria

(4) Tube 2nd set of labels

(5) Laser radiation label for collimator AL04 II-D

Only for deliveries to China

or


(6) Laser radiation label for collimator RFU


(7) Color coding label

(only when WRCC or Remote Interface is configured)

### 7.1.3 Labels on the wall stand

(1) Wall stand ID label

(2) Unit revision label

This product complies with DHHS regulations  
21 CFR Subchapter J applicable at date of  
manufacture.  
Manufactured August 2018  
Bairong Shanghai Medical Equipment Ltd  
278 Zhongyuan Road, 201310 Shanghai  
China

## (3) DHHS label

### SECOND SET OF LABELS DETECTOR

## (4) Detector 2nd set of label



## (5) Core static ID label



## (6) 25 kg label

## (7) Type B label

## (8) Iontomat label (only the right label)



## (9) Wireless detector position labels

For Core XL, only one label is pasted on the table (left) and wall stand (right).


(10) Anti-grip label

### 7.1.4 Labels on the patient table

(1) Wireless detector position labels

For Core XL, only one label is pasted on the table (left) and wall stand (right).

(2) Iontomat label for table with wireless detector


(3) Warning label of crush leg


(4) Anti-grip label

(5) Patient table ID label


(6) Unit revision label

(7) Iontomat label for table with fixed detector

### 7.1.5 Labels on the generator ON/OFF console

#### WARNING!

This X-ray unit may be dangerous to patient and operator unless safe exposure factors, operating instructions and maintenance schedules are observed.


(1) Radiation emission warning label for deliveries to USA

(2) DHHS label for deliveries to USA

### 7.1.6 Labels on the generator cover

(1) Caution label

(2) Refer to manual label

(3) High voltage label

**SECOND SET OF LABELS  
HV-GENERATOR**

(4) Generator 2nd set of label

**SECOND SET OF LABELS  
PSU**

(5) PSU 2nd set of label

## 7.2 Technical data

All technical data represent typical values unless specific tolerances are given.

Siemens Healthineers reserves the right to change the design and specifications without notice.

### 7.2.1 System

|                                                                                     |                        |                                                                                              |                                                                                                                                                                 |
|-------------------------------------------------------------------------------------|------------------------|----------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Installation data                                                                   |                        | The entire system is powered via a three phase voltage connection.                           |                                                                                                                                                                 |
| Power connection                                                                    |                        | 3-phase, 380 V/400 V/440 V (50/60 Hz), 480 V (60 Hz) (± 10%)                                 |                                                                                                                                                                 |
| Power consumption                                                                   |                        | Max. 127 kVA (80 kW)                                                                         |                                                                                                                                                                 |
| Standby                                                                             |                        | 1.2 kVA                                                                                      |                                                                                                                                                                 |
|                                                                                     |                        |                                                                                              |                                                                                                                                                                 |
| Internal line resistance                                                            | Max. R <sub>Line</sub> |                                                                                              |                                                                                                                                                                 |
| U <sub>Line</sub>                                                                   | 55 kW                  | 65 kW                                                                                        | 80 kW                                                                                                                                                           |
| 380 V                                                                               | 0.15 Ω                 | 0.15 Ω                                                                                       | 0.10 Ω                                                                                                                                                          |
| 400 V                                                                               | 0.17 Ω                 | 0.17 Ω                                                                                       | 0.11 Ω                                                                                                                                                          |
| 440 V                                                                               | 0.20 Ω                 | 0.20 Ω                                                                                       | 0.14 Ω                                                                                                                                                          |
| 480 V                                                                               | 0.24 Ω                 | 0.24 Ω                                                                                       | 0.16 Ω                                                                                                                                                          |
|                                                                                     |                        |                                                                                              |                                                                                                                                                                 |
| Ambient conditions (operation):<br>Examination room                                 | Temperature range      | +10°C to +35°C                                                                               |                                                                                                                                                                 |
|                                                                                     | Relative humidity      | 20% to 75%, non-condensing                                                                   |                                                                                                                                                                 |
|                                                                                     | Barometric pressure    | 700 hPa to 1060 hPa                                                                          |                                                                                                                                                                 |
|                                                                                     | Height                 | ≤ 3000 m above sea level                                                                     |                                                                                                                                                                 |
| Ambient conditions (operation):<br>Imaging system                                   | Temperature range      | 0°C to +35°C                                                                                 |                                                                                                                                                                 |
|                                                                                     | Relative humidity      | 20% to 75%, non-condensing                                                                   |                                                                                                                                                                 |
| Ambient conditions (Transport<br>and storage)                                       | Temperature range      | -20°C to +70°C                                                                               |                                                                                                                                                                 |
|                                                                                     | Relative humidity      | 10% to 95%, non-condensing                                                                   |                                                                                                                                                                 |
|                                                                                     | Barometric pressure    | 700 hPa to 1060 hPa                                                                          |                                                                                                                                                                 |
| Radio interference suppression/EMC                                                  |                        | IEC 60601-1-2                                                                                |                                                                                                                                                                 |
| Protection class                                                                    |                        | The system belongs to protection class I with a type B applied part according to IEC 60601-1 |                                                                                                                                                                 |
|  |                        | Protection against ingress of water:                                                         | <ul style="list-style-type: none"><li>• Foot switches: IPx8</li><li>• MAX wI-D: IP43 or higher</li><li>• Core XL: IPx3</li><li>• Rest of system: IPx0</li></ul> |

|                                                                                 |                                                                                                                                                                           |                                               |
|---------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------|-----------------------------------------------|
|                                                                                 | Attention: To avoid the risk of electrical shock, a protective conductor must be implemented.                                                                             |                                               |
| Explosion protection<br>(AP, APG after IEC 60601-1-3:2013;<br>IEC 60601-1:2012) | The system is not intended for operation in explosion-endangered areas (for example highly flammable mixtures of anaesthesia gases with air or oxygen or nitrous oxide).  |                                               |
| Operating mode                                                                  | Continuous operating with intermittent loading                                                                                                                            |                                               |
| Weight                                                                          | Bucky wall stand                                                                                                                                                          | Approx. 240 kg                                |
|                                                                                 | X-ray tube support                                                                                                                                                        | With 3 m rail: approx. 300 kg                 |
|                                                                                 |                                                                                                                                                                           | With 4 m rail: approx. 320 kg                 |
|                                                                                 | Imaging station                                                                                                                                                           | With normal screen: approx. 14.5 kg           |
|                                                                                 |                                                                                                                                                                           | With touch screen (optional): approx. 15.4 kg |
|                                                                                 | Generator cabinet                                                                                                                                                         | Approx. 190 kg                                |
|                                                                                 | Patient table                                                                                                                                                             | Approx. 320 kg                                |
| Degree of contamination                                                         | Class 2, according to IEC 60601-1-3:2013; IEC 60601-1:2012                                                                                                                |                                               |
| Total inherent filtration                                                       | 2.5 mm Al at 70 kV/2.5 mm Al HVL (DIN EN 61267, RQR5) (combined with collimator Al equivalent 1.0 mm Al and tube Al equivalent 1.5 mm Al) + DAP Al equivalent < 0.5 mm Al |                                               |
| DAP chamber                                                                     | Dose area product is measured by DAP meter at the collimator                                                                                                              |                                               |
|                                                                                 | Tolerance of DAP values                                                                                                                                                   | ± 35%                                         |
| AEC chambers                                                                    | AEC chambers attenuation equivalent                                                                                                                                       | 0.7 mm Al                                     |
|                                                                                 | Tolerance of AEC values                                                                                                                                                   | ± 3%                                          |
| Applied parts                                                                   | <ul style="list-style-type: none"> <li>• Table</li> <li>• Compression device</li> <li>• Wall stand</li> <li>• Portable detector</li> </ul>                                |                                               |

## 7.2.2 Components

### Patient table

|                        |                                                                                                                                            |
|------------------------|--------------------------------------------------------------------------------------------------------------------------------------------|
| Table height           | 51.5 cm to 90.0 cm; total lift 38.5 cm (tabletop)                                                                                          |
| Tabletop dimensions    | <ul style="list-style-type: none"> <li>• Standard tabletop: 233 cm x 80 cm</li> <li>• Short tabletop (optional): 213 cm x 80 cm</li> </ul> |
| Tabletop Al equivalent | ≤ 0.7 mm Al (with 100 kV / 3.6 mm Al HVL; IEC 60601-2-54)                                                                                  |
| Patient weight         | Max. load 300 kg                                                                                                                           |

|                              |                                                                                                                                                                                                                                                                                                               |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Longitudinal travel          | <ul style="list-style-type: none"> <li>• Standard tabletop: <math>\pm 44</math> cm</li> <li>• Short tabletop (optional): <math>\pm 34</math> cm</li> </ul>                                                                                                                                                    |
| Transverse travel            | $\pm 14$ cm ( $\pm 5$ mm tolerance)                                                                                                                                                                                                                                                                           |
| Detector cover range         | $\geq 100$ cm                                                                                                                                                                                                                                                                                                 |
| Max. patient coverage        | About 190 cm with standard tabletop                                                                                                                                                                                                                                                                           |
| Grid                         | Stationary, Pb 13/92, $f_0 = 115$ cm                                                                                                                                                                                                                                                                          |
|                              | Stationary, Pb 13/40; $f_0 = 115$ cm                                                                                                                                                                                                                                                                          |
| Tabletop - detector distance | $\leq 73$ mm                                                                                                                                                                                                                                                                                                  |
| Maximum radiation field      | <ul style="list-style-type: none"> <li>• Maximum symmetrical radiation field at a distance of 1150 mm from the focus,</li> <li>• With table tray for MAX wi-D: 426 mm x 355 mm</li> <li>• With table tray for Core XL: 426 mm x 426 mm</li> <li>• With table tray for Core static: 426 mm x 426 mm</li> </ul> |

### Bucky wall stand

|                                      |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Travel range (central beam to floor) | From 31.5 cm to 175 cm                                                                                                                                                                                                                                                                                                                                                                                                      |
| Detector unit                        | Tiltable from $-20^\circ$ to $+90^\circ$ with detent at $0^\circ$ and $+90^\circ$                                                                                                                                                                                                                                                                                                                                           |
| Detector cover - detector distance   | <ul style="list-style-type: none"> <li>• With MAX wi-D: <math>\leq 4.2</math> cm</li> <li>• With Core XL: <math>\leq 4.2</math> cm</li> <li>• With Core static: <math>\leq 4.5</math> cm</li> <li>• With a cassette: <math>\leq 4.6</math> cm</li> </ul>                                                                                                                                                                    |
| X-ray absorption                     | $\leq 0.6$ mm Al (at 100 kV/3.6 mm Al HVL; IEC 60601-2-54)                                                                                                                                                                                                                                                                                                                                                                  |
| Anti scatter grid                    | <ul style="list-style-type: none"> <li>• Stationary grid, Pb 13/92, <math>f_0 = 140</math> cm</li> <li>• Stationary grid, Pb 13/40, <math>f_0 = 140</math> cm</li> </ul> Optional stationary grids: <ul style="list-style-type: none"> <li>• Pb 13/92, <math>f_0 = 115</math> cm; <math>f_0 = 180</math> cm; <math>f_0 = 300</math> cm</li> <li>• Pb 13/40, <math>f_0 = 115</math> cm; <math>f_0 = 180</math> cm</li> </ul> |

### X-ray tube support

|                         |                                        |        |
|-------------------------|----------------------------------------|--------|
| Horizontal travel range | Longitudinal                           | 352 cm |
|                         | Transverse with 3 m trolley:           | 215 cm |
|                         | Transverse with 4 m trolley:           | 350 cm |
| Vertical travel range   | 180 cm $\pm$ 2 cm, manual or motorized |        |

|                                           |                                                                                              |
|-------------------------------------------|----------------------------------------------------------------------------------------------|
| Minimum focus-ceiling distance            | 83 cm                                                                                        |
| Rotation range around the vertical axis   | -154° to +182°, manual<br>Detents at -90°; 0°; +90°; +180°                                   |
| Rotation range around the horizontal axis | ± 137° ± 2°, motorized tilting;<br>± 120° ± 2°, manual tilting;<br>Detents at -90°; 0°; +90° |

### X-ray tube

|                                                                                |                                                   |
|--------------------------------------------------------------------------------|---------------------------------------------------|
| Maximum exposure voltage (IEC 60613)                                           | 150 kV                                            |
| Focal spot nominal value (IEC 60336)                                           | • 0.6<br>• 1.2                                    |
| Radiographic anode input power (IEC 60613)                                     | 34 kW/80 kW                                       |
| Nominal anode input power (IEC 60613) (thermal anode reference output = 300 W) | 22 kW/54 kW (50/60Hz)<br>33 kW/78 kW (150/180Hz)  |
| Nominal anode input power (thermal anode reference output = 0 W)               | 29 kW/73 kW (50/60Hz)<br>45 kW/107 kW (150/180Hz) |
| Anode angle                                                                    | 12°                                               |
| Maximum heat storage capacity of tube housing                                  | 1,000 kJ (1,350 kHU)                              |
| Total filtration (IEC 60601-1-3)                                               | ≥ 2.5 mm Al/75 kV                                 |
| Leakage radiation (IEC 60601-1-3)                                              | (at 150kV at 1 m distance) ≤ 0.8 mGy/h (450 W)    |

### Collimator

|                            |                                                                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Inherent filtration        | • Collimator AL04 II-D: 1mm Al at 70 kV<br>• Collimator RFU: 1mm Al at 75 kV                                                                                                                                                  |
| Full field light localizer | Very efficient high power LED technology, high energy efficiency enabling low-noise design without external cooling system; long lifetime approx. 100,000 h; timer functionality; laser line light localizer (can be covered) |
| Copper prefilter           | Without filter; 0.1 mm; 0.2 mm; 0.3 mm; motorized                                                                                                                                                                             |
| Rotation                   | ±45°                                                                                                                                                                                                                          |
| Collimator control         | Manual or motorized                                                                                                                                                                                                           |

## Detectors

| MAX wi-D                                | pixium 3543 EZh                                                                                                                                                                                       |
|-----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detector technology                     | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology                                                                                                                    |
| Dimensions (active area)                | 34.8 cm x 42.4 cm                                                                                                                                                                                     |
| Dimensions with detector housing        | 44 cm x 46.1 cm x 1.9 cm                                                                                                                                                                              |
| Active detector matrix                  | 2350 x 2866                                                                                                                                                                                           |
| Pixel size                              | 148 $\mu$ m                                                                                                                                                                                           |
| Semiconductor material                  | Amorphous silicon (a-Si)                                                                                                                                                                              |
| Scintillator                            | Cesium iodide (CsI)                                                                                                                                                                                   |
| Acquisition depth                       | 16-bit                                                                                                                                                                                                |
| Data transmission                       | < 2 s preview; < 4 s full image                                                                                                                                                                       |
| DQE in %; 2 $\mu$ Gy (RQA5) (IEC 62220) | 70% at 0.05 lp/mm<br>51% at 1 lp/mm<br>42% at 2 lp/mm<br>29% at 3 lp/mm<br>19% at Nyquist                                                                                                             |
| MTF in % (RQA5) (IEC 62220)             | 63% at 1 lp/mm<br>35% at 2 lp/mm<br>19% at 3 lp/mm<br>12% at Nyquist                                                                                                                                  |
| Thickness                               | 19 mm                                                                                                                                                                                                 |
| Weight                                  | 3.3 kg                                                                                                                                                                                                |
| Max. load capacity                      | <ul style="list-style-type: none"> <li>300 kg with patient recumbent on it (300 kg is only for functional availability rather than full performance)</li> <li>100 kg with patient standing</li> </ul> |
| Battery                                 | Lithium-ion, rechargeable, exchangeable                                                                                                                                                               |
| Charging time                           | 3 h                                                                                                                                                                                                   |
| Battery operation time                  | <ul style="list-style-type: none"> <li>Up to 1,050 images</li> <li>6.5 hours in usual operation</li> </ul>                                                                                            |
| Charging location                       | Table detector tray, wall stand detector tray, battery charger (optional)                                                                                                                             |
| Clip-on grid                            | <ul style="list-style-type: none"> <li>Pb 5/85, <math>f_0 = 115</math> cm</li> <li>Pb 15/80, <math>f_0 = 115</math> cm</li> </ul>                                                                     |

|                                          |                                                                                                                                 |
|------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| MAX wi-D                                 | pixium 3543 EZh                                                                                                                 |
| WLAN Standard                            | IEEE 802.11 b/g/n, 2 x 2 mimo, WPA2/AES Encryption, EAP/TLS support                                                             |
| Core XL                                  | Mars1717VS                                                                                                                      |
| Detector technology                      | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology                                              |
| Dimensions (active area)                 | 42.6 cm x 42.6 cm                                                                                                               |
| Dimensions with detector housing         | 46.1 cm x 46.1 cm x 1.57 cm                                                                                                     |
| Active detector matrix                   | 3070 x 3070                                                                                                                     |
| Pixel size                               | 139 $\mu$ m                                                                                                                     |
| Semiconductor material                   | Amorphous silicon (a-Si)                                                                                                        |
| Scintillator                             | Cesium iodide (CsI)                                                                                                             |
| Acquisition depth                        | 16-bit                                                                                                                          |
| Data transmission                        | < 3 s preview; < 7 s full image                                                                                                 |
| DQE in %; 2 $\mu$ Gy (RQA5), (IEC 62220) | 80% at 0.05 lp/mm<br>65% at 1 lp/mm<br>53% at 2 lp/mm<br>34% at 3 lp/mm<br>21% at Nyquist                                       |
| MTF in % (RQA5), (IEC 62220)             | 64% at 1 lp/mm<br>34% at 2 lp/mm<br>18% at 3 lp/mm<br>13% at Nyquist                                                            |
| Thickness                                | 15.7 mm                                                                                                                         |
| Weight                                   | 4.2 kg                                                                                                                          |
| Max. load capacity                       | <ul style="list-style-type: none"> <li>• 150 kg with patient recumbent on it</li> <li>• 100 kg with patient standing</li> </ul> |
| Battery                                  | Lithium-ion, rechargeable, exchangeable                                                                                         |
| Charging time                            | 4 h                                                                                                                             |
| Battery operation time                   | <ul style="list-style-type: none"> <li>• Up to 950 images</li> <li>• 7.5 hours in usual operation</li> </ul>                    |
| Charging location                        | Table detector tray, wall stand detector tray, battery charger (optional)                                                       |
| Clip-on grid                             | Pb 5/85, $f_0 = 115$ cm                                                                                                         |

|                                          |                                                                                           |
|------------------------------------------|-------------------------------------------------------------------------------------------|
| <b>Core XL</b>                           | <b>Mars1717VS</b>                                                                         |
| WLAN Standard                            | IEEE 802.11 b/g/n, 2 x 2 mimo, WPA2/AES Encryption, EAP/TLS support                       |
| <b>Core static</b>                       | <b>Venu1717X</b>                                                                          |
| Detector technology                      | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology        |
| Dimensions (active area)                 | 42.6 cm x 42.6 cm                                                                         |
| Dimensions with detector housing         | 46.0 cm x 46.0 cm x 1.55 cm                                                               |
| Active detector matrix                   | 3070 x 3070                                                                               |
| Pixel size                               | 139 $\mu$ m                                                                               |
| Semiconductor material                   | Amorphous silicon (a-Si)                                                                  |
| Scintillator                             | Cesium iodide (CsI)                                                                       |
| Acquisition depth                        | 16-bit                                                                                    |
| Data transmission                        | < 3 s preview; < 5 s full image                                                           |
| DQE in %; 2 $\mu$ Gy (RQA5), (IEC 62220) | 80% at 0.05 lp/mm<br>65% at 1 lp/mm<br>53% at 2 lp/mm<br>34% at 3 lp/mm<br>21% at Nyquist |
| MTF in % (RQA5), (IEC 62220)             | 64% at 1 lp/mm<br>34% at 2 lp/mm<br>18% at 3 lp/mm<br>13% at Nyquist                      |

## Generator

All technical data prove typical values, if specific tolerances are not indicated.

### Polydoros RFX

According to IEC 60601-2-7 and GB 9706.3

|                                |                   |
|--------------------------------|-------------------|
| High voltage form              | Multipulse        |
| Generator frequency            | $\geq 100$ kHz    |
| Exposure voltage               | 40 kV - 150 kV    |
| Maximum current                | 50 A              |
| Maximum kW's value in AEC mode | 56 kW's $\pm 5$ % |

**Polydoros RFX**


|                               |                                                                                                                                                                                                                                                           |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| mAs range                     | <ul style="list-style-type: none"> <li>• 0.5 mAs to 800 mAs for 55 kW/65 kW</li> <li>• 0.5 mAs to 1000 mAs for 80 kW</li> </ul>                                                                                                                           |
| Minimum exposure time         | <ul style="list-style-type: none"> <li>• 1-point technique: 1 ms at 100 kV</li> <li>• 2-point technique: 2 ms at 100 kV</li> <li>• 3-point technique: 20 ms at 100 kV</li> </ul>                                                                          |
| Tolerance                     | <ul style="list-style-type: none"> <li>• kV tolerance <math>\pm 5\%</math></li> <li>• mAs tolerance not greater than <math>\pm (10\% + 0.2 \text{ mAs})</math></li> <li>• ms tolerance not greater than <math>\pm (10\% + 1 \text{ ms})</math></li> </ul> |
| Nominal electric power        | <ul style="list-style-type: none"> <li>• 55 kW : 100 kV, 55 mAs, 100 ms (550 mA)</li> <li>• 65 kW : 100 kV, 65 mAs, 100 ms (650 mA)</li> <li>• 80 kW : 100 kV, 80 mAs, 100 ms (800 mA)</li> </ul>                                                         |
| Max. electric power           | <ul style="list-style-type: none"> <li>• 55 kW : 100 kV, 55 mAs, (550 mA)</li> <li>• 65 kW : 100 kV, 65 mAs, (650 mA)</li> <li>• 80 kW : 100 kV, 80 mAs, (800 mA)</li> </ul>                                                                              |
| Max. tube current and voltage | <ul style="list-style-type: none"> <li>• 55 kW: 68 kV at 800 mA</li> <li>• 65 kW: 81 kV at 800 mA</li> <li>• 80 kW: 80 kV at 1,000 mA</li> </ul>                                                                                                          |
| mA value at nominal 150 kV    | <ul style="list-style-type: none"> <li>• 55 kW: 350 mA</li> <li>• 65 kW: 400 mA</li> <li>• 80 kW: 509 mA</li> </ul>                                                                                                                                       |

**Image Station**

|                             |                         |                                                                                                                                 |
|-----------------------------|-------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| Computer system             | Computer                | Intel® Core™ i5-6500TE Processor, 6 M Cache, up to 3.30 GHz, 1 x DDR4 8 GB memory, 2 x 500 G SATA HDD, 4 x USB 2.0, 2 x USB 3.0 |
|                             | Operating system        | Windows® 10 LTSB 64 bit or above                                                                                                |
| Archiving and documentation | System storage capacity | 10,000 images                                                                                                                   |
|                             | Background functions    | Preprocessing and data transfer is done in the background.                                                                      |
|                             | Archiving               | USB, DVD or CD-R read/write in DICOM 3.0 format (or packed format)                                                              |

|  |                      |                                                                                                                                                     |
|--|----------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------|
|  | CD-R storage:        | <ul style="list-style-type: none"> <li>• DICOM 3.0 format</li> <li>• Packed format</li> <li>• DICOM send</li> <li>• DICOM query/retrieve</li> </ul> |
|  | Export               | In JPG, BMP, or TIFF format via USB/DVD/CD-R                                                                                                        |
|  | Documentation        | <ul style="list-style-type: none"> <li>• DICOM Print</li> <li>• Virtual Film Sheet</li> </ul>                                                       |
|  | DICOM Worklist, MPPS | Data exchange via HIS/RIS                                                                                                                           |
|  | Paper printing       | Printing of lists and images to a network paper printer                                                                                             |
|  |                      |                                                                                                                                                     |

### 7.2.3 Network information

- (1) #3 NIC  
 (2) #2 NIC  
 (3) Imaging system (Windows 10 LTSB 2016)  
 (4) #1 NIC  
 (5) DICOM TCP Port 104  
 (6) Switch on customer site  
 (7) Any RIS or HIS or PACS

MULTIX Impact C uses the Windows firewall to restrict access.

The following ports need to be opened in the hospital network firewall:

| Application              | Functionality | Network - MULTIX Impact C (In-going) | Network - MULTIX Impact C (Out-going) |
|--------------------------|---------------|--------------------------------------|---------------------------------------|
| Windows operating system | DNS           | None                                 | 53 (UDP, TCP)                         |
|                          | DHCP          | None                                 | 67 (UDP)                              |
|                          | NTP           | None                                 | 123 (UDP)                             |

| Application     | Functionality      | Network - MULTIX Impact C (In-going) | Network - MULTIX Impact C (Out-going) |
|-----------------|--------------------|--------------------------------------|---------------------------------------|
| DICOM           | DICOM              | 104 (TCP)                            | Configurable by customer* (TCP)       |
| Remote Services | MNP                | 8226, 13001 (TCP)                    | 8227, 8228, 12061 (TCP)               |
|                 | Teamviewer         | 10080, 10081 (TCP)                   | None                                  |
|                 | Transfermanager    | 49152-65534 (TCP)                    | 20-21 (TCP)                           |
|                 | Remote Web Service | 443 (TCP)                            | None                                  |

\*Ports for DICOM communication are opened for configured peer nodes, e.g. PACS or RIS or HIS.

### 7.3 Notice concerning electromagnetic compatibility (EMC)

Medical electrical equipment needs special precautions regarding EMC. EMC information provided in the accompanying documents must be followed where appropriate.

#### WARNING

The portable detector generates an electromagnetic field so high that it disturbs life-supporting device.

#### **Interference with life supporting device resulting in possible malfunction**

- ◆ Keep a safety separation distance of more than 30 cm between the portable detector and the life-supporting device.
- ◆ In case of interference with other equipment increase the distance between the interfering devices.

#### CAUTION

Radio interference from the system affects nearby equipment

#### **Risk of malfunction of nearby equipment.**

- ◆ Reorient or relocate the equipment or system or shield its location.
- ◆ Test the operation of the equipment in its new location.

The system is suitable for professional healthcare facility environment according to standard IEC 60601-1-2:2014.

Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30 cm (12 inches) to any part of the system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

Fixed product cabling, which can not be removed by the users, is not listed. This cabling is part of the product and was present at all EMC-measurements. Without this cabling there is no complete functionality of the product.

Deviations or additions to this document are added to the product specific documents.

The use of accessories, transducers and cables other than those specified, with the exception of transducers and cables sold by the manufacturer of the product as replacement parts for internal components, may result in increased emission or decreased immunity of the product.

The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required) this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measures, such as relocating or reorienting the equipment.

The system should not be used adjacent to or stacked with other equipment and that if adjacent or stacked use is necessary, the equipment or product should be observed to verify normal operation in the configuration in which it will be used.

Precautions to take if the use location is near (e.g. less than 1.5 km from) AM, FM or TV broadcast antennas.

The essential performance of this equipment is defined in particular medical electrical equipment standard, such as IEC 60601-2-43, IEC 60601-2-54.

Medical electrical equipment is subject to special precautions regarding EMC. These systems must be installed and put into service according to the EMC information provided in the accompanying documents.

### 7.3.1 Guidance and manufacturer's declaration - electromagnetic emissions

The system is intended for use in the electromagnetic environment specified below. The customer or the user of the product should assure that it is used in such an environment.

| Emission test                                            | Compliance     | Electromagnetic environment - guidance                                                                                                                                                                                                                                                                                                                                                          |
|----------------------------------------------------------|----------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| RF emissions<br>CISPR 11                                 | Group 1        | The system uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.<br><br>The system is suitable for use in all establishments other than domestic and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes. |
| RF emissions<br>CISPR 11                                 | Class A        |                                                                                                                                                                                                                                                                                                                                                                                                 |
| Harmonic emissions<br>IEC 61000-3-2                      | Not applicable |                                                                                                                                                                                                                                                                                                                                                                                                 |
| Voltage fluctuations/ flicker emissions<br>IEC 61000-3-3 | Not applicable |                                                                                                                                                                                                                                                                                                                                                                                                 |

### 7.3.2 Guidelines and manufacturer's declaration - electromagnetic immunity

The system is intended for use in an electromagnetic environment as specified below. The customer or the user of the system should ensure that it is operated in such an environment.

| Immunity test                                                                                          | IEC 60601 Test level                                                                                                      | Compliance level                                                          | Electromagnetic environment - Guidelines                                                                                                                                                                                                                            |
|--------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------|---------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Electrostatic discharge (ESD)<br>IEC 61000-4-2                                                         | $\pm 8$ kV contact<br>$\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air                                                 | $\pm 8$ kV contact<br>$\pm 2$ kV, $\pm 4$ kV, $\pm 8$ kV, $\pm 15$ kV air | Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.                                                                                                                       |
| Electrical fast transient/ burst<br>IEC 61000-4-4                                                      | $\pm 2$ kV for power port<br>$\pm 1$ kV for input/ output port<br>100 kHz repetition frequency                            | The same as IEC 60601 test level                                          | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                 |
| Surge<br>IEC 61000-4-5                                                                                 | $\pm 0.5$ kV, $\pm 1$ kV line to line<br>$\pm 0.5$ kV, $\pm 1$ , $\pm 2$ kV line to ground                                | The same as IEC 60601 test level                                          | Mains power quality should be that of a typical commercial or hospital environment.                                                                                                                                                                                 |
| Voltage dips, short interruptions and voltage variations on power supply input lines<br>IEC 61000-4-11 | 0% $U_T$ ; 0.5 cycle<br>(0°, 45°, 90°, 135°, 180°, 225°, 270° and 315°)<br>0% $U_T$ ; 1 cycle<br>70% $U_T$ ; 25/30 cycles | Not applicable                                                            | The mains power quality should be that of a typical commercial or hospital environment. If the user of the system requires continued operation during power mains interruptions, it is recommended that the system be powered from an uninterruptible power supply. |
|                                                                                                        | 0% $U_T$ ; 250/300 cycles                                                                                                 | < 5% $U_T$ for 0.5 sec (> 95% dip in $U_T$ )                              |                                                                                                                                                                                                                                                                     |
| Power frequency magnetic field<br>IEC 61000-4-8                                                        | 30 A/m<br>50 Hz or 60 Hz                                                                                                  | 30 A/m                                                                    | Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.                                                                                                                           |

NOTE:  $U_T$  is the a.c. mains voltage prior to application of the test level.

| Immunity test | IEC 60601 Test level | Compliance level | Electromagnetic environment - Guidelines                                                                                                                                                                                                                           |
|---------------|----------------------|------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|               |                      |                  | Portable and mobile RF communications equipment should be used no closer to any part of the product than the recommended separation distance calculated from the equation applicable to the frequency of the transmitters.<br><br>Recommended separation distance: |

| Immunity test                                                                              | IEC 60601 Test level                                                                        | Compliance level                           | Electromagnetic environment - Guidelines                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------|--------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Conducted RF<br>IEC 61000-4-6                                                              | 3 V<br>150 kHz-80 MHz<br>6 V; in ISM bands between<br>150 kHz and 80 MHz<br>80% AM at 1 kHz | 3 V<br>6 V; in ISM bands                   | $d=1.2 \sqrt{P}$                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Radiated RF<br>IEC 61000-4-3<br>Proximity fields from RF wireless communications equipment | 3 V/m<br>80 MHz to 2.7 GHz<br>80% AM at 1 kHz<br>See 8.10 of<br>IEC 60601-1-2-2014          | 3 V/m<br>See 8.10 of<br>IEC 60601-1-2-2014 | $d=1.2 \sqrt{P}$ for 80 MHz to 800 MHz<br>$d=2.3 \sqrt{P}$ for 800 MHz to 2.7 GHz<br><p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m).</p> <p>Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey<sup>1</sup>, should be less than the compliance level in each frequency range<sup>2</sup>.</p> <p>Interference may occur in the vicinity of equipment marked with the following</p> <p>symbol: </p> |

Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

<sup>1</sup> Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To be able to assess the electromagnetic environment with regard to the stationary transmitters, an electromagnetic site study should be considered. If the measured field strength at the location where the system is used exceeds the applicable RF compliance level specified above, the system should be observed to verify normal operation. Should unusual performance features be observed, additional measures (such as change in orientation or change of site of the system) may be necessary.

<sup>2</sup> Above the frequency range of 150 kHz to 80 MHz the field strength should be less than 3 V/m.

### 7.3.3 Recommended separation distance between portable and mobile RF communications equipment and the system

The system is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the system can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the system as recommended below, according to the maximum output power of the communications equipment.

| Nominal power of the transmitter [W] | Separation distance according to frequency of transmitter [m] |                                       |                                        |
|--------------------------------------|---------------------------------------------------------------|---------------------------------------|----------------------------------------|
|                                      | 150 kHz to 80 MHz<br>$d=1.2 \sqrt{P}$                         | 80 MHz to 800 MHz<br>$d=1.2 \sqrt{P}$ | 800 MHz to 2.7 GHz<br>$d=2.3 \sqrt{P}$ |
| 0.01                                 | 0.12                                                          | 0.12                                  | 0.23                                   |
| 0.1                                  | 0.38                                                          | 0.38                                  | 0.73                                   |
| 1                                    | 1.2                                                           | 1.2                                   | 2.3                                    |
| 10                                   | 3.8                                                           | 3.8                                   | 7.3                                    |
| 100                                  | 12                                                            | 12                                    | 23                                     |

For transmitters rated at a maximum output power not listed above, the recommended separation distance  $d$  in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where  $P$  is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Remark 1: At 80 MHz and 800 MHz, the higher frequency range applies.

Remark 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

### 7.3.4 RF transmitters

| Name                                   | Frequency range                                          | Bandwidth                                                                          | Modulation type                                                                                                     | Effective radiated power                                                               |
|----------------------------------------|----------------------------------------------------------|------------------------------------------------------------------------------------|---------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------|
| Wireless Remote Control Console (WRCC) | 2402 - 2480 MHz/2 MHz step                               | 2402 - 2480 MHz/2 MHz step                                                         | GFSK                                                                                                                | 27 dB<br>31.6 mA (Max)                                                                 |
| Wireless detector MAX wi-D             | 2400 - 2483.5 MHz,<br>5150 - 5350 MHz,<br>5725- 5850 MHz | 2400 - 2483.5 MHz,<br>5150 - 5350 MHz,<br>5725- 5850 MHz                           | BPSK/QPSK/16QAM/<br>64QAM/256QAM/<br>DBPSK/DQPSK/CCK                                                                | 33 dBm (Max)<br>20 dBm (EIRP) (Max)                                                    |
| Wireless detector Core XL              | 2400 - 2483.5 MHz,<br>5150 - 5250 MHz,<br>5725- 5850 MHz | 2.4GHz: 20MHz, 40MHz<br>5.1GHz: 20MHz, 40MHz, 80MHz<br>5.8GHz: 20MHz, 40MHz, 80MHz | 11b: CCK, DQPSK, DBPSK, DSSS<br>11n: QPSK, BPSK, 16QAM, 64QAM, OFDM<br>11ac: QPSK, BPSK, 16QAM, 64QAM, OFDM, 256QAM | 802.11b: 16dBm<br>802.11a: 13dBm<br>802.11g: 14dBm<br>802.11n: 13dBm<br>802.11ac: 8dBm |
| iPad mini A2133                        | 2400-2483.5 MHz,<br>5150-5350 MHz,<br>5725-5850 MHz      | 80 MHz, 20 MHz, 8 MHz, 3 MHz                                                       | BPSK/QPSK/16QAM/<br>256QAM/DBPSK/<br>DQPSK/CCK GFSK<br>π/4-DQPSK 8DPSK                                              | 33 dBm (Max)<br>23 dBm (EIRP) (Max)<br>20 dBm (EIRP) (Max)                             |

## 8 Glossary

|       |                                                |
|-------|------------------------------------------------|
| ACSS  | Automatic Collimation Size Sensing             |
| AEC   | Automatic Exposure Control                     |
| BWS   | Bucky Wall stand                               |
| DICOM | Digital Imaging and Communications in Medicine |
| DMG   | Dead Man's Grip                                |
| EMC   | Electromagnetic compatibility                  |
| MPPS  | Modality Performed Procedure Step              |
| OGP   | Organ Program                                  |
| RHA   | Rotation about the horizontal axis             |
| ROI   | Region of Interest                             |
| RVA   | Rotation about the vertical axis               |
| SID   | Source Image Distance                          |
| TOD   | Table-Object Distance                          |
| TUI   | Touch User Interface                           |
| UPS   | Uninterruptible Power Supply                   |
| WRCC  | Wireless Remote Control Console                |

**1, 2, 3 ...**

- 2D camera
  - adjusting ROI for Ortho 153

**A**

- Abnormal movements 27
- Accessories 157
  - description 159
  - orientation 159
  - safety information 158
- Accessory rails
  - multileaf collimator 69
- acquisition
  - Ortho 155
- ACSS 137
- ACSS = Automatic Collimation Size Sensing 137
- AEC 133
- AEC = Automatic Exposure Control 133
- Application
  - capabilities 53
  - patient table 72
  - wall stand 75
- Attaching the detector 127
- Automatic exposure control
  - checking 91
- Auxiliary devices 157

**B**

- BABIX holder 170, 178
  - positioning 180
  - usage 170
- Barcode reader 216
- Battery
  - changing 122, 124
- Battery charger 212, 214
  - status LEDs 213
  - troubleshooting 213
- Brake knob 72
- Bucky wall stand
  - Ortho 149

**C**

- Calibration
  - detector 128
- camera
  - adjusting ROI for Ortho 153
- Cassette holder
  - attaching 194
  - danger points 36
  - lateral 192
- Cassette holder without retainer 194
- Ceiling-mounted tube support 58
- Centering
  - checking 89
- Centering display
  - TUI 65

- Charging
  - wireless detector 125
- Cleaning 162
  - safety notes 47
- Clinical benefits 8
- Clip-on grids 198
  - attaching 198
  - removing 199
- Collimation 135
  - Automatic Collimation Size Sensing 137
  - checking 89
  - restoring 138
  - TUI 65
- Collimator mode
  - TUI 65
- Collision protection 27
- Combination
  - with other products/components 46
- Compensation filters 200
  - attachment 202
- Compression Belt 184
  - cleaning and disinfection 185
  - description 186
  - location of labels 190
  - operation 187
  - safety check 185
  - safety information 184, 185
  - standards 191
  - technical data 189
- Conditions of use 9
- Contraindications 8
- Control field
  - TUI 61
- Control room
  - overview 55
- Copyright 20
- Core XL 77
- Covers 41
- Cu filter 68, 134

**D**

- Daily checks 88
- Danger of crushing 27
- Danger points 33
- Danger zones 33
- DAP
  - display test 92
- Data security 15
- Deadman function 27
- Detaching the detector 128
- Detector 77
  - handling 117
  - management 126
- Detector calibration 128
- Detector holder 181
  - attaching 196
  - danger points 36

- Detector holder, lateral 196
- Detector sharing 128
- Detector tracking
  - TUI 63
- Detector tray 73
  - danger points 35
  - positioning 104
- Detector unit
  - movements 109
- DICOM conformity 20
- Disinfection 162
  - safety notes 47
- Displays
  - safety information 39

**E**

- Electric shock
  - protection 41
- EMC = Electromagnetic compatibility 42
- Emergency operation 129
- Emergency procedures 32
- Emergency SHUTDOWN button 31
- Emergency STOP
  - patient table 74
- Emergency STOP buttons 29
  - locations 29
- Equipment
  - cleaning 49
  - damage 44
  - safety measures 41
- Examination
  - checks 90
  - safety information 37
- Examination room
  - overview 54
- Explosion protection 16
- exposure
  - Ortho 155
- Exposure
  - free exposure 144
- Exposure release
  - with hand switch 57
- Exposures
  - onto the patient table 138
  - onto the wall stand 142

**F**

- FD calibration 128
- FD management 126
- FD protector 199
- Filter holder 204
  - operation 204
- Fire
  - safety 16
- Floating tabletop 27
- Foot kick switches 74
- Foot switch 211

Free exposure 144  
 Front panel markings 77  
 Functional check 88

## G

Generator  
   damage 45  
   ON/OFF console 56  
   TUI 63

Grid 105  
   removing 113  
   wall holder 200

Grid status  
   TUI 65

Grip  
   locations 36

## H

Hand grip  
   attaching 167  
   locations 36

Hand grip strip 164  
   attaching 165

Hand grip, lateral 166

Hand switch  
   exposure release 57

Hardcopy cameras  
   safety information 39

Horizontal projection  
   onto the mobile detector 146  
   onto the wall stand 142

## I

Illustrations 13

Image orientation 37

Indications for use 8

Information  
   Operator Manual 11

Installation 26

Intended purpose 8

Intended use 8

Intended users 9

Intercom system  
   in control room 217

Interfaces 46

Introduction 8

IONTOMAT  
   automatic exposure control 133

IT-security 20

## L

Labels  
   location 218

Laser line light localizer 70

Lateral cassette holder  
   danger points 36

Lateral detector holder  
   danger points 36

Laws 14

LCD  
   UPS 84

Light localizer 70, 135  
   checking 89  
   safety information 39  
   switching on 136

Localizer 70

Location  
   of labels 218

## M

Maintenance  
   periodic 21  
   plan 21  
   preventive 23  
   safety inspection 22

Mattress 163  
   attaching 164

MAX wi-D 77

Maximum load 43

Maximum weight 43

Measuring tape 69

Mobile detector holder CR 195

Modifications 26

Monthly checks 91

movements  
   buttons and indications 60

Movements  
   checking 89  
   detector unit 109  
   patient table 101

Multileaf collimator 67  
   accessory rails 69  
   rotating 99

Multipurpose / Ortho Stand  
   technical data 182

## N

Names 12

Network security 17

## O

Oblique projection  
   onto the patient table 141  
   onto the wall stand 143

OEM 19

OEM = Original Equipment Manufacturer 19

OFF 81

ON 81

ON/OFF console 56

Operation 161  
   on battery power 86

Operator Manual  
   information 11

organ program  
   for Ortho 150

Ortho 148  
   adjusting ROI with camera 153  
   adjusting ROI with start and end position 152  
   on patient table 149  
   on wall stand 149  
   organ program 150  
   positioning on the Ortho stand 151  
   positioning patient on patient table 151  
   preparing 150  
   releasing exposures 155

Ortho Stand 171

Other products/components 46

Overhead patient hand grip 168

## P

Paper roll holder 164

Parameters 12

Patient  
   positioning at the wall stand 132  
   positioning on patient table 151  
   positioning on the table 131  
   visual contact 38

Patient group 9

Patient information  
   TUI 61

Patient orientation 37

Patient positioning  
   safety information 37

Patient positioning mattress 163

Patient size  
   TUI 62

Patient support 199

Patient table 72  
   danger points 34  
   Emergency STOP 74  
   exposures 138  
   foot kick switches 74  
   height adjustment 101  
   movements 101

Pediatric Use 147

Periodic maintenance 21

Physical functionality 9

Pictograms 12

Positioning  
   patient on patient table 151  
   the patient at the wall stand 132  
   the patient on the table 131

Potential equalization 41

Power failure 32

Power outlet 41

Power supply 41

Preventive maintenance 23

Product  
   description 53

Product description 160

Protection  
 against electric shock 41  
 Protection class 41  
 Purpose  
 of the product 8

## R

Radiation protection 40  
 Red emergency STOP buttons 29  
 locations 29  
 Regulations 14  
 Remote software update 20  
 Repair 26  
 Replacement  
 of safety-relevant parts 25  
 Requirements  
 for diagnosis 38  
 Residual risks 9  
 RHA 95  
 RHA = rotation about the horizontal axis 95  
 Risk of crushing 161  
 ROI  
 for Ortho with camera 153  
 for Ortho with start and end  
 position 152  
 Room lighting 38  
 Ruler 177  
 RVA 95  
 RVA = rotation about the vertical axis 95

## S

Safety  
 general 14  
 personal 27  
 Safety check 88, 163  
 Safety inspection 22  
 Safety-relevant parts 25, 43  
 Scattered radiation grid 105  
 Service life  
 product 43  
 SID display 99  
 SID tracking 98  
 TUI 63  
 SID/SOD measurement 69  
 Side effects 9  
 Smart virtual Ortho 153  
 Smart Virtual Ortho  
 SmartOrtho 148  
 SmartOrtho 152  
 Software  
 copyright 20  
 remote update 20  
 safety information 17  
 third-party 19  
 STOP buttons 29  
 locations 29

Swivel stool 176  
 System configurations 54  
 System overview  
 control room 55  
 examination room 54

System status  
 TUI 65

## T

Table height  
 adjusting 101  
 Tabletop  
 damage 44  
 positioning 102  
 Technical data 228  
 Multipurpose / Ortho Stand 182  
 Technical description 163  
 Technical documents 26  
 Test images 38  
 Third-party software 19  
 Three-field template  
 storage 205  
 usage 205  
 Three-field template, set 205  
 TOD  
 for Ortho 151  
 Top alignment  
 TUI 63  
 Tracking  
 SID 98  
 tube 97  
 Training 14  
 Trolley UM mobile patient table 159  
 Tube position  
 TUI 65  
 Tube support 58  
 Tube tracking 97  
 TUI 63  
 Tube unit  
 directions for moving and rotating 93  
 positioning 94  
 TUI 59  
 displaying help 61  
 layout 59  
 movements 60  
 TUI = Touch User Interface 59

## U

Unit movements  
 safety information 27  
 Units  
 opening 41  
 UPS 83  
 description 83  
 LCD 84  
 operation 85  
 troubleshooting 87  
 UPS = Uninterruptible Power Supply 83

Usability 10  
 Use  
 of the product 8  
 User training 14

## V

Value statements 13  
 Vertical projection  
 onto the patient table 139  
 onto the pulled-out detector 140  
 onto the wall stand 143  
 Visual contact 38

## W

Wall charger 215  
 Wall holder for grids 200  
 Wall stand 75  
 danger points 35  
 exposures 142  
 horizontal projection 142  
 oblique projection 143  
 operating elements 75  
 safety information 108  
 vertical projection 143  
 Warning signs 45  
 Warning sound 45  
 Wireless detector 77  
 battery change 122, 124  
 charging 125  
 description 77  
 detector tray 73  
 handling 117  
 maintenance 120  
 management 126  
 messages 79  
 recommendations 116  
 safety information 115  
 sharing 128  
 WLAN 120  
 WLAN  
 operation 120  
 Workplace selection  
 TUI 62  
 WRCC 206  
 battery exchange 211  
 charging 210  
 displays 208  
 marking signs 211  
 operating elements 209  
 operation 207  
 technical data 211  
 WRCC = Wireless Remote Control Console 206

## X

X-ray tube support  
 danger points 34

Original language: English

---

**Siemens Healthineers Headquarters**  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany  
Phone: +49 9131 84-0  
[siemens.com/healthineers](https://www.siemens.com/healthineers)

**Manufacturer**  
Siemens Shanghai Medical Equipment  
Ltd.  
278 Zhou Zhu Road  
201318 Shanghai  
P.R.China

**EU Authorized Representative**  
Siemens Healthcare GmbH  
Henkestr. 127  
91052 Erlangen  
Germany

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The screenshot shows a login interface with two input fields: 'User Name' and 'Password'. The 'User Name' field has a dropdown arrow on the right. Below the fields are two buttons: 'Log On' and 'Emergency'.

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The screenshot shows the 'User Name' dropdown menu open, displaying a list of users. The 'Password' field and the 'Log On' and 'Emergency' buttons are also visible.

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The screenshot shows the 'User Name' field populated with 'Mgr' and the 'Password' field filled with dots. The 'Log On' and 'Emergency' buttons are also visible.

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| Exam Protocol |    |                        |    |            |                   |           |                 |         |                     |                          | MULTIX Impact |      |               |       |           |
|---------------|----|------------------------|----|------------|-------------------|-----------|-----------------|---------|---------------------|--------------------------|---------------|------|---------------|-------|-----------|
| Patient Name  |    | HER_201811120812411*** |    |            | Request ID        |           |                 |         |                     | Number of images         |               | 1    |               |       |           |
| Date of Birth |    |                        |    |            | Accession No.     |           | 201811120812411 |         |                     |                          |               |      |               |       |           |
| Patient ID    |    | 201811120812411        |    |            | Study Description |           | Emergency       |         |                     | Dose Area Prod.(mGy*cm2) |               | 0.00 |               |       |           |
| Patient Sex   |    | D                      |    |            | Ref. Physician    |           |                 |         |                     | Ref. Air Kerma(mGy)      |               | 0    |               |       |           |
| Ward          |    |                        |    |            | Institution Name  |           |                 |         |                     |                          |               |      |               |       |           |
| Isolation no. | kV | ma                     | mm | Dist. Lev. | MField            | Cu.Filter | Format          | Density | Date and Time       | Application              | DAP           | EN   | Ref.Air Kerma | Focus | Grid      |
| 1             | 40 | 5                      | 0  |            | M                 | NONE      |                 | 0       | 2018-11-12 10:44:07 |                          | 0.00          | 0.57 | 0             | SMALL | Universal |

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Operator Manual – Imaging System  
VA20

# Legend

|                                                                                             |                                                                                                                                                                                                                                                                                                                                                                                                                                               |
|---------------------------------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|            | Indicates a hint<br>Is used to provide information on how to avoid operating errors or information emphasizing important details                                                                                                                                                                                                                                                                                                              |
|            | Indicates the solution of a problem<br>Is used to provide troubleshooting information or answers to frequently asked questions                                                                                                                                                                                                                                                                                                                |
|            | Indicates a list item                                                                                                                                                                                                                                                                                                                                                                                                                         |
|            | Indicates a prerequisite<br>Is used for a condition that has to be fulfilled before starting a particular operation                                                                                                                                                                                                                                                                                                                           |
|            | Indicates a one-step operation                                                                                                                                                                                                                                                                                                                                                                                                                |
|            | Indicates steps within operating sequences                                                                                                                                                                                                                                                                                                                                                                                                    |
| <i>Italic</i>                                                                               | Is used for references and for table or figure titles                                                                                                                                                                                                                                                                                                                                                                                         |
|           | Is used to identify a link to related information as well as previous or next steps                                                                                                                                                                                                                                                                                                                                                           |
| <b>Bold</b>                                                                                 | Is used for screen outputs of the system                                                                                                                                                                                                                                                                                                                                                                                                      |
| Menu > Menu Item                                                                            | Is used for the navigation to a certain submenu entry                                                                                                                                                                                                                                                                                                                                                                                         |
| Text                                                                                        | Is used to identify inputs you need to provide                                                                                                                                                                                                                                                                                                                                                                                                |
| Orange                                                                                      | Is used to emphasize particularly important sections of the text                                                                                                                                                                                                                                                                                                                                                                              |
|  CAUTION | CAUTION<br>Used with the safety alert symbol, indicates a hazardous situation which, if not avoided, could result in minor or moderate injury or material damage.<br>CAUTION consists of the following elements: <ul style="list-style-type: none"> <li>• Information about the nature of a hazardous situation</li> <li>• Consequences of not avoiding a hazardous situation</li> <li>• Methods of avoiding a hazardous situation</li> </ul> |
|  WARNING | WARNING<br>Indicates a hazardous situation which, if not avoided, could result in death or serious injury.<br>WARNING consists of the following elements: <ul style="list-style-type: none"> <li>• Information about the nature of a hazardous situation</li> <li>• Consequences of not avoiding a hazardous situation</li> <li>• Methods of avoiding a hazardous situation</li> </ul>                                                        |

|          |                                                     |           |
|----------|-----------------------------------------------------|-----------|
| <b>1</b> | <b>Safety notes</b>                                 | <b>6</b>  |
| 1.1      | Usability                                           | 6         |
| 1.2      | Federal law                                         | 6         |
| 1.3      | Scope of applicability                              | 6         |
| 1.4      | General safety notes                                | 6         |
| 1.5      | Image failures                                      | 7         |
| 1.6      | Remote software update                              | 8         |
| 1.7      | Language of the user interface                      | 9         |
| 1.8      | Information about this document                     | 9         |
| 1.8.1    | Operator Manual                                     | 9         |
| 1.8.2    | Options                                             | 9         |
| 1.8.3    | Illustrations                                       | 10        |
| 1.8.4    | Value statements                                    | 10        |
| 1.8.5    | Names in images                                     | 10        |
| <b>2</b> | <b>Welcome to the MULTIX Impact C Image Station</b> | <b>11</b> |
| 2.1      | Application                                         | 11        |
| 2.2      | Acquisition mode                                    | 11        |
| <b>3</b> | <b>Startup and shutdown</b>                         | <b>12</b> |
| 3.1      | Switching on and starting                           | 12        |
| 3.2      | Unlock the disk                                     | 13        |
| 3.3      | Login and logoff                                    | 13        |
| 3.3.1    | Logging in                                          | 14        |
| 3.3.2    | Logging off                                         | 16        |
| 3.4      | Switching off the system                            | 16        |
| <b>4</b> | <b>Operation of the imaging system</b>              | <b>17</b> |
| 4.1      | Screen layout                                       | 17        |
| 4.2      | Online help                                         | 18        |
| 4.2.1    | Accessing the Operator Manual                       | 18        |
| 4.2.2    | Accessing the Online Help                           | 19        |
| 4.3      | Status information                                  | 19        |
| 4.3.1    | General status                                      | 20        |
| 4.3.2    | Detector status                                     | 21        |
| 4.3.3    | Status in the patient list                          | 22        |
| 4.3.4    | System messages                                     | 23        |
| 4.3.5    | Software update available                           | 24        |
| 4.3.6    | Status of image transfers                           | 24        |
| 4.4      | Study lists and patient registration                | 26        |
| 4.4.1    | Study lists                                         | 26        |
| 4.4.2    | Accessing the study lists                           | 26        |
| 4.4.3    | Searching for studies/patients                      | 28        |
| 4.4.4    | Selecting a study/patient                           | 30        |
| 4.4.5    | Patient registration                                | 31        |
| 4.4.6    | Correcting study data                               | 34        |
| 4.4.7    | Deleting studies                                    | 35        |
| 4.4.8    | Protecting/unprotecting studies                     | 36        |
| 4.4.9    | FD management                                       | 36        |
| 4.4.10   | Quality control study                               | 39        |
| 4.5      | Image acquisition                                   | 41        |
| 4.5.1    | General information                                 | 41        |
| 4.5.2    | Information on acquisition using a detector         | 44        |
| 4.5.3    | Checking/changing organ programs                    | 45        |

|        |                                                     |     |
|--------|-----------------------------------------------------|-----|
| 4.5.4  | Acquisition guidance                                | 47  |
| 4.5.5  | Checking/changing acquisition parameters            | 48  |
| 4.5.6  | Organ program and acquisition parameters on the TUI | 54  |
| 4.5.7  | Releasing an exposure                               | 55  |
| 4.5.8  | Exposure index                                      | 56  |
| 4.5.9  | Accepting or rejecting images                       | 58  |
| 4.5.10 | Repeating an exposure                               | 60  |
| 4.5.11 | Image recovery                                      | 60  |
| 4.5.12 | Closing a study                                     | 61  |
| 4.5.13 | MPPS close                                          | 63  |
| 4.5.14 | Ortho operation                                     | 64  |
| 4.6    | Stitch mode                                         | 71  |
| 4.6.1  | Viewing an Ortho image                              | 71  |
| 4.6.2  | Composing                                           | 73  |
| 4.7    | Image postprocessing                                | 76  |
| 4.7.1  | Processing images                                   | 76  |
| 4.7.2  | Adding a laterality label                           | 78  |
| 4.7.3  | Adding texts                                        | 78  |
| 4.7.4  | Changing the text size                              | 78  |
| 4.7.5  | Deleting text items                                 | 79  |
| 4.7.6  | Windowing principle                                 | 79  |
| 4.7.7  | Windowing an image                                  | 80  |
| 4.7.8  | Inverting the gray scale of an image                | 81  |
| 4.7.9  | Cropping an image                                   | 81  |
| 4.7.10 | Zooming an image                                    | 82  |
| 4.7.11 | Panning an image                                    | 83  |
| 4.7.12 | Flipping an image                                   | 84  |
| 4.7.13 | Rotating an image                                   | 84  |
| 4.7.14 | Setting the position and laterality information     | 85  |
| 4.7.15 | Undoing changes                                     | 87  |
| 4.7.16 | Saving changes                                      | 87  |
| 4.7.17 | Setting an image flavor                             | 88  |
| 4.7.18 | Customized flavor                                   | 90  |
| 4.7.19 | Using the Rad tool                                  | 92  |
| 4.8    | Image evaluation                                    | 94  |
| 4.8.1  | Viewing images                                      | 95  |
| 4.8.2  | Measuring angles                                    | 96  |
| 4.8.3  | Measuring distances                                 | 96  |
| 4.8.4  | Deleting measurements                               | 97  |
| 4.8.5  | Calibrating images                                  | 97  |
| 4.8.6  | Drawing arrows                                      | 100 |
| 4.8.7  | Storing an image                                    | 101 |
| 4.8.8  | Marking images                                      | 101 |
| 4.8.9  | Displaying the exam protocol                        | 102 |
| 4.8.10 | Sending/archiving images                            | 102 |
| 4.9    | Image hardcopy                                      | 103 |
| 4.9.1  | Arranging film sheets                               | 104 |
| 4.9.2  | Changing the image size on the film sheet           | 106 |
| 4.9.3  | Changing the image view on the film sheet           | 108 |
| 4.9.4  | Printing images                                     | 108 |
| 4.10   | Study/patient administration                        | 109 |
| 4.10.1 | Displaying the images of studies                    | 109 |
| 4.10.2 | Importing studies/images                            | 110 |
| 4.10.3 | Exporting studies/images                            | 112 |

|   |          |                                             |            |
|---|----------|---------------------------------------------|------------|
|   | 4.10.4   | Sending/archiving studies/images            | 114        |
|   | 4.10.5   | Retrieving studies/images                   | 116        |
|   | 4.10.6   | Managing studies in the Recycle Bin         | 116        |
| ■ | <b>5</b> | <b>Administration of the imaging system</b> | <b>118</b> |
|   | 5.1      | User settings                               | 118        |
|   | 5.1.1    | Site information                            | 119        |
|   | 5.1.2    | General settings                            | 121        |
|   | 5.1.3    | Image processing properties                 | 122        |
|   | 5.1.4    | Examination properties                      | 124        |
|   | 5.1.5    | Image display properties                    | 125        |
|   | 5.1.6    | Patient registration configuration          | 127        |
|   | 5.1.7    | DICOM properties                            | 128        |
|   | 5.1.8    | Disk management                             | 131        |
|   | 5.1.9    | Settings for Ortho                          | 132        |
|   | 5.2      | User management                             | 134        |
|   | 5.2.1    | User roles and privileges                   | 134        |
|   | 5.2.2    | Changing the password                       | 134        |
|   | 5.2.3    | Setting up users, roles, and permissions    | 135        |
|   | 5.3      | Organ program and procedure editor          | 140        |
|   | 5.3.1    | General information                         | 140        |
|   | 5.3.2    | Opening the OGP Editor                      | 142        |
|   | 5.3.3    | The OGP Editor                              | 143        |
|   | 5.3.4    | Editing procedures                          | 143        |
|   | 5.3.5    | Procedure data                              | 144        |
|   | 5.3.6    | Editing organ programs                      | 146        |
|   | 5.3.7    | Organ program parameters                    | 149        |
|   | 5.3.8    | Configuration parameters                    | 155        |
|   | 5.3.9    | Closing the OGP Editor                      | 157        |
|   | 5.4      | Service functions                           | 158        |
|   | 5.4.1    | Clearing study and image data               | 158        |
|   | 5.4.2    | Detector calibration                        | 159        |
|   | 5.4.3    | FD reconnect                                | 163        |
|   | 5.4.4    | RIS mapping                                 | 163        |
|   | 5.4.5    | Statistics                                  | 168        |
|   | 5.4.6    | Audit trail                                 | 169        |
|   | 5.4.7    | Test images                                 | 170        |
|   | 5.4.8    | Tube load                                   | 172        |
|   | 5.4.9    | Providing an image for service              | 172        |
|   | 5.4.10   | Providing log files for service             | 172        |
|   | 5.4.11   | Service login                               | 173        |
|   | 5.4.12   | Remote assistance                           | 173        |
|   | 5.4.13   | Remote service access                       | 174        |
|   | 5.4.14   | DICOM nodes                                 | 175        |
|   | 5.4.15   | FD manager                                  | 180        |
|   | 5.4.16   | Sensitivity calibration                     | 180        |
|   | 5.4.17   | QR codes                                    | 181        |
|   | 5.4.18   | HIPAA                                       | 181        |
| ■ | <b>6</b> | <b>Glossary</b>                             | <b>183</b> |
| ■ |          | <b>Index</b>                                | <b>184</b> |

# 1 Safety notes

Welcome to the MULTIX Impact C Operator Manual.

Please pay attention to the relevant safety information.

For general notes on safety see the Safety part of the Operator Manual of the MULTIX Impact C.

## 1.1 Usability

### BEFORE USING THE SYSTEM

This system should be used only by employees adequately trained in the use of this equipment.

Before using this system, the operator should be thoroughly acquainted with the instructions for use and safety recommendations provided in this manual.

Failure to follow the instructions for use and safety recommendations provided in this manual can cause serious injury to the patient, to the operator or to other persons.

## 1.2 Federal law

21 CFR 801.109 (b)(1)

Caution: Federal law restricts this device to sale by or on the order of a physician.

## 1.3 Scope of applicability

This document applies to MULTIX Impact C with software version VA20 and higher.

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**Risk of malfunctions posing a hazard to patients and equipment**

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## 1.5 Image failures

- ◆ Always pay attention to the image quality.

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**CAUTION**

Radiation has been released and no image appears on the monitor of the apparently ready imaging system.

**Risk of radiation without a stored image**

- ◆ If a black image appears (image text changes but image is underexposed), make a test image before continuing patient examinations.

If no image appears (no change on the display):

- Do not release any further X-rays.
- If the portable detector was used, attempt image recovery.  
If image recovery is unsuccessful, go to the next step.
- Shut down the system and call Siemens Healthineers Service.

**CAUTION**

An acquisition has been released and no image appears on the monitor of the apparently ready imaging system.

**Risk of X-ray radiation without a stored image**

- ◆ Do not release any further acquisitions.
- ◆ Shut down the system and call Siemens Healthineers Service.

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**Risk of unnecessary radiation exposure**

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## 1.6 Remote software update

When remote software update is enabled, software packages are installed from Siemens Healthineers:

- Microsoft "hotfixes"
- System software updates

When a software package is available, you are prompted to accept the installation.

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## 1.7 Language of the user interface

Your system can be configured by the Siemens Healthineers Service in one of the following languages (user interface and keyboard layout):

- English
- German
- French
- Spanish
- Chinese
- Japanese

Some system messages are in English only.

## 1.8 Information about this document

### 1.8.1 Operator Manual

The Operator Manual for your system consists of the System Operator Manual, this manual, and the Addenda.

Please read the Addenda to the Operator Manual.

### 1.8.2 Options

This document describes all system features of all system variants.

The complete system is described with all options and components that have been released. Possible options have not been marked specially.

Particular options or components may not be available for specific systems.

The quotation text of your order is the sole reference for the functional scope of your system.

If your system does not have a specific feature, please contact your local sales representative.

### 1.8.3 Illustrations

All illustrations of equipment and of the program user interface shown in this document are examples only. The functions actually available depend on the type of system, the installed options, and the current configuration.

Other differences in detail may occur in your system due to constant development and improvement of the system.

Reproduction of images can cause loss of detail. Pictures in this document do not therefore provide any indication of image quality.

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Values shown in pictures of the software user interface have no clinical meaning.

Please only set the values preset in the organ programs provided or the values recommended by experienced application specialists.

### 1.8.5 Names in images

All names of patients shown in images are purely fictional. Any similarities with existing persons are entirely coincidental.

## 2 Welcome to the MULTIX Impact C Image Station

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With only a basic knowledge of the operation of PCs and Windows®, you can quickly learn how to operate the imaging system.

The main purpose of the digital imaging system is to manage patients (studies), acquire, display, postprocess, film, and archive X-ray images.

### 2.1 Application

The imaging system is suitable for image acquisition with the MULTIX Impact C system.

It can be used for many types of radiological examinations.

Please also refer to the Operator Manual of the MULTIX Impact C.

### 2.2 Acquisition mode

The imaging system provides the Digital Radiography (RAD) acquisition mode using a solid-state flat detector (FD) or cassette.

(→ Page 49 *Workplaces and exposure techniques*)

## 3 Startup and shutdown

The imaging system is switched on and off together with the whole system.

- (→ Page 12 *Switching on and starting*)
- (→ Page 13 *Login and logoff*)
- (→ Page 16 *Switching off the system*)

Do not try to switch the system on again immediately after switch-off. Wait until the imaging system is shut down (dark monitor). Then, wait a few seconds before switching on again.

### 3.1 Switching on and starting

- 1 Press the "ON" button.

The system is powered up.

- 2 Wait until the application program is displayed on the monitor.

- The system automatically runs a self-test when it is switched on.
- The operating system of the computer is started.
- The application program is started.

If the display of logon banner is enabled, the logon banner is displayed.

- 3 Click OK.

The logon banner text can be configured.

(→ Page 120 *Login banner*)

The login screen is displayed.

The system is now ready to use.

In case of a fault

- ◆ Shut down the system and turn it on again.
- If the fault remains:
- ◆ Call Siemens Healthineers Service.

Power failure

- ◆ Switch off the entire system and then switch it on again. The last image might be lost.
- If a fault is indicated when the system starts up again:
- ◆ Call Siemens Healthineers Service.

### 3.2 Unlock the disk

The hard disk containing the patient data and image can be encrypted by the Siemens Healthineers service.

✓ The data encryption is activated.

1 Enter the encryption password.

2 Click **Unlock**.

The login screen is displayed.

A recovery password is saved during the data encryption by the Siemens Healthineers service. Keep it in a safe place, e.g. a USB stick.

If the encryption password is lost, forgotten, or unavailable, you can use the recovery password to unlock the disk.

You can change the encryption password after a successful recovery.

Do not backup the system when the hard disk has been encrypted.

### 3.3 Login and logoff

HIPAA provides user management assigning each user different rights. You must log in using your login name and password. This prevents the system, patient data and images from misuse by unauthorized persons.

Your system administrator must create user accounts and assign groups and privileges to them.

Depending on the rights given to a user by the administrator, not all patients can be seen and some functions of the imaging system may not be available.

### 3.3.1 Logging in

#### Local login

- 1 Enter your login name.

– or –

Select a user account from the list.

- 2 Enter your password.
  - 3 Click Log On.
- The **Pre-registered Patients** list is displayed.  
(→ Page 134 *Changing the password*)

#### Domain login

A registered domain user can log into the imaging system if an active directory is available to the system.

- ✓ Domain login has been turned on.

The screenshot shows a login interface with a dark background. At the top right, there are two buttons: 'Domain' and 'Local'. Below these are three input fields labeled 'Domain name', 'User Name', and 'Password'. At the bottom right, there are two buttons: 'Log On' and 'Emergency'.

- 1 Enter your domain name.
- 2 Enter the registered user name and the password.
- 3 Click **Log On**.

The **Pre-registered Patients** list is displayed.

The default user role for a domain user is **Tech**.

After first login, the administrator can change the user role for the domain user in the **Service** mode.

(→ Page 134 *User management*)

However, the password of the domain user can only be changed in the domain server.

Local login is still possible by turning the switch to **Local**.

(→ Page 14 *Local login*)

### Performing an emergency login

If you forget your password or you do not have an account for this system and must perform an emergency examination, you can perform an emergency login.

- ◆ Click **Emergency** on the login screen.

A new study for the emergency patient is registered. The name of the emergency patient is set to: **ER\_**current date and time. The new study is passed to the **Examination** mode. You need to select the proper organ program to continue the study.

When logged in as an emergency user, it is only possible to perform acquisition. Only emergency patients can be seen.

Patients acquired in emergency mode are assigned by the system to no group, i.e. everybody can see them.

### 3.3.2 Logging off

If another user takes over or you want to leave the room you should log off for safety reasons.

- ◆ Click the user icon on the top right.

The login window is displayed.

If the system is not operated for a configurable time period, it will be locked automatically via the screen saver.

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- 1 Terminate the current examination or review.

- 2 Press the "OFF" button.

When shutting down the system during an active image transfer, a message about the incomplete transfer job will occur.

- 3 You select **Finish Transfer**.

The imaging system monitors the status of the transfer jobs, and when the last job is completed, the system is shut down. This may take some time.

– or –

You select **Immediate Shutdown**.

The transfer jobs are aborted. The transfer jobs are kept in the list and can be restarted when the system is started up again. The system is shut down.

The application program is terminated.

The operating system of the computer shuts down (dark monitor).

The entire system including all connected devices is switched off.

#### Power off

Do not switch off the system using the main switch before the operating system of the computer has been shut down.

- ◆ Wait until the imaging system has been shut down (dark monitor).

Only when the monitor is dark, you may switch off using the main switch.

## 4 Operation of the imaging system

- (→ Page 17 *Screen layout*)
- (→ Page 26 *Study lists and patient registration*)
- (→ Page 41 *Image acquisition*)
- (→ Page 76 *Image postprocessing*)
- (→ Page 94 *Image evaluation*)
- (→ Page 103 *Image hardcopy*)
- (→ Page 109 *Study/patient administration*)


The imaging system functions can be operated from the control room and some functions also from the examination room.

### 4.1 Screen layout

Screen layout in Patient mode

- (1) Access bar
- (2) Tabs of the patient list
- (3) Display area for patient lists, images and dialogs
- (4) Status and system messages bar
- (5) Time
- (6) Available functions

Screen layout in Examination, Postprocessing, Stitch, View, and Print mode

- (1) Access bar
- (2) Tabs of the task cards for mode selection
- (3) Information area
- (4) Status and system messages bar
- (5) Time
- (6) Available functions
- (7) Display area for images and instructions

If the system is not operated for a configurable time period, the screen content will be hidden via the screen saver.

(→ Page 121 *General settings*)

## 4.2 Online help


You can display the Operator Manual in the Knowledge Gateway or you can enable context-sensitive Online Help.

### 4.2.1 Accessing the Operator Manual

- 1 Click the ? icon top right.

The Knowledge Gateway is displayed accordingly.

#### Example

- 2 Select the required language of the Operator Manual and double-click the file name. The Operator Manual is displayed.

For more information about the functions in the Knowledge Gateway, please click the **? Help** icon in the Knowledge Gateway.

### 4.2.2 Accessing the Online Help

- ◆ Move the cursor on a button and press the **F1** key. The related chapter in the Online Help is displayed. If the cursor is not positioned on a button, the Knowledge Gateway is displayed.

The Online Help can be configured in one of the following languages:

- English
- German
- French
- Spanish
- Chinese

### 4.3 Status information

Status information is given on the screen of the imaging system via icons in the status bar bottom left and text information bottom right.

Example of status information right: current time

Example of status icons bottom left

- (→ Page 20 *General status*)
- (→ Page 21 *Detector status*)
- (→ Page 22 *Status in the patient list*)
- (→ Page 23 *System messages*)
- (→ Page 24 *Status of image transfers*)
- (→ Page 174 *Remote service access*)

### 4.3.1 General status

| Icon                                                                                            | Short description                                                                                            |
|-------------------------------------------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------|
|  disconnected  | Network status                                                                                               |
|  connected     |                                                                                                              |
|  ≤90%          | Free hard disk space<br>The disk fill level for blocking acquisition is configurable. #DROC Disk management# |
|  ≤87.5%      |                                                                                                              |
|  ≤75%        |                                                                                                              |
|  ≤50%        |                                                                                                              |
|  ≤25%        |                                                                                                              |
|  ≤10%        |                                                                                                              |
|  TUBE        | Tube heat units                                                                                              |
|  n           | (→ Page 23 <i>System messages</i> )                                                                          |
|              | Software update is available.                                                                                |
|  standby     | (→ Page 174 <i>Remote service access</i> )                                                                   |
|  ongoing     |                                                                                                              |
|  in progress | DICOM send job status                                                                                        |
|  failed      |                                                                                                              |

| Icon                                                                                          | Short description |
|-----------------------------------------------------------------------------------------------|-------------------|
|  in progress | Print job status  |
|  failed      |                   |

When you hold the mouse cursor over an icon, status information will be displayed.

### 4.3.2 Detector status

| Icon                                                                                                               | Short description                                                                                                            |
|--------------------------------------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------------------|
|  0...40%                          | Wireless signal strength of the selected detector                                                                            |
|  40...60%                         |                                                                                                                              |
|  60...80%                        |                                                                                                                              |
|  80...100%                      |                                                                                                                              |
|  within calibration temperature | Detector temperature                                                                                                         |
|  out of calibration temperature |                                                                                                                              |
|  <-10°C                         | Detector temperature is not within required range.<br>Acquisitions are blocked.                                              |
|  >46°C                          |                                                                                                                              |
|  0...10%                        | Battery charging condition of the selected detector<br><br>Acquisitions are blocked when the battery level is less than 10%. |
|  10...20%                       |                                                                                                                              |
|  20...40%                       |                                                                                                                              |
|  40...60%                       |                                                                                                                              |
|  60...80%                       |                                                                                                                              |
|  80...100%                      |                                                                                                                              |

| Icon                                                                              | Short description                 |
|-----------------------------------------------------------------------------------|-----------------------------------|
|  | Detector calibration has expired. |

Pay attention to the status icons.

The detector specific icons in the status area always refer to the currently selected detector.

When the battery icon is yellow, detector battery is low. Charge the detector battery or replace it by a fully charged battery.

### 4.3.3 Status in the patient list

In the patient list, status icons indicate the current study status:

| Icon                                                                                | Short description                                                                                                |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------------------------------------------------|
|   | The study is scheduled.<br>The study comes in via HIS/RIS and its status has not been switched to a local study. |
|  | The study is in progress.                                                                                        |
|  | The study is completed and closed.                                                                               |
|  | The patient information has been changed.                                                                        |
|  | The study has been imported.                                                                                     |
|  | The study has been exported.                                                                                     |
|  | The study has been retrieved from network.                                                                       |
|  | The study has been burned into CD/DVD.                                                                           |
|  | The study is protected from modification and deletion.                                                           |
|  | The study is not protected from modification and deletion.                                                       |
|  | Archiving is in progress.                                                                                        |
|  | Archiving succeeded.                                                                                             |


| Icon                                                                              | Short description          |
|-----------------------------------------------------------------------------------|----------------------------|
|  | Archiving failed.          |
|  | DICOM send is in progress. |
|  | DICOM send succeeded.      |
|  | DICOM send failed.         |
|  | Print is in progress.      |
|  | Print succeeded.           |
|  | Print has failed.          |

#### 4.3.4 System messages

You can check system messages at any time, e.g. in case of an error.

The message type is indicated by a colored dot:

- Green dot: general message
- Yellow dot: warning message
- Red dot: error message
- ◆ Click this icon in the status bar.

(Example message)

### Deleting messages

When a message is no more relevant, you can delete it:

- ◆ Click this icon.

— or —

Click **Clear Warning**.

All warning messages are cleared in the list.

### Refreshing messages

You can refresh an error message after the remedy has been performed.

- ◆ Click this icon.

#### 4.3.5 Software update available

When a software update is available for your imaging system, a confirmation dialog appears after the system is switched on.

- 1 Click **Install** to start the installation.

The system is blocked during software update.

— or —

Click **Defer** to start the installation at a later time.

An icon indicates in the status bar:

You can click the icon in the status bar once you are ready to start the installation.

A confirmation dialog appears after the successful installation.

- 2 Click **OK**.

The system will reboot automatically.


#### 4.3.6 Status of image transfers

All image transfers run in the background.

### Checking the job status

- 1 Click this icon.

- 2 Click **Send** or **Print**, if necessary.

### Starting a job once again

You can restart an image transfer which has not been completed:

- 1 Select the images in the appropriate transfer list.  
You can use the **Ctrl** or **Shift** key to select multiple images.
- 2 Click **Retry**.

### Deleting a job

Please check that the images have been successfully transferred before you delete them.

You can abort image transfers in certain phases during preparation of the job and between jobs:

- 1 Select the images in the appropriate transfer list.  
You can use the **Ctrl** or **Shift** key to select multiple images.
- 2 Click **Delete**.

### Updating the job status list

- ◆ Click **Refresh** to update the job status list.

## 4.4 Study lists and patient registration

All patient data with the personal data of a patient, the acquisitions and other examination and image data is stored in studies on the hard disk of the imaging system.

Patient data is stored in compliance with the DICOM standard.

Patient data is arranged in study lists in **Patient mode**.

(→ Page 26 *Study lists*)

Before starting an examination, a patient registration must be performed, either via HIS/RIS or manually.

(→ Page 31 *Patient registration*)

### 4.4.1 Study lists

You can access patient data via the following study lists in **Patient mode**:

- You will enter the personal data of the patient and create a new study in the **Pre-registered Patients** list (work list). Alternatively, patient data (a study or a procedure step) may come in via HIS/RIS (hospital/radiological information system).
- All completed and closed studies (images, series) are stored under the selected study entry in the **Examined Patients** list.
- When an examination is finished you will archive the study data via network, CD/DVD/USB memory, or hardcopy. You may then delete the image data to regain space for new patients.

Studies retrieved from a DICOM node are stored both in the **Examined Patients** and **Archiving Information** list.

- If activated, the **Recycle Bin** list contains deleted studies. These can be restored at any time.

### Single or multi study mode

The study lists can be configured for **Single study mode** or **Multi study mode**.

(→ Page 121 *General settings*)

- The **Single study mode** treats each study individually.  
By clicking on a study in **Single study mode** just one study is selected.
- The **Multi study mode** treats the studies of a patient together.  
By clicking a study in **Multi study mode**, all studies of the patient with the same date are selected.

### 4.4.2 Accessing the study lists

#### Selecting a study list

- 1 Click the **Patient** tab, if necessary.

- 2 Click the tab of the required study list.

You can configure the default position of the mouse cursor after entering Patient mode.

(→ Page 121 *General settings*)

If the tabs of Archiving Information and Recycle Bin are not visible, ask your administrator to assign the permissions.

(→ Page 136 *Assigning permissions to roles*)

## Changing the displayed study entries

In **Examined Patients** and **Recycle Bin** tab, you can assign the study entry number displayed on one page.

- ◆ Enter an desired patients number in the entry of **Patient Per Page** and press the **Enter** key.

Click the left or right arrow to show previous or next page.

## Scrolling a study list

- ◆ Scroll using the scroll bar or the mouse wheel.

## Changing the sorting sequence of study entries

You can change the sorting sequence by the column of **Name**, **ID**, **Sex**, **Age**, **Description**, and **Date/Time**.

- ◆ Click in the header of a column you require.

The study list is then sorted alphanumerically by this column.

The study entries are sorted by the patients' **First Name** if you click in the header of **Name**.

The new sorting sequence only applies to the study entries on current displayed page.

You can configure the displayed items, position, and sizes in the study lists.  
(→ Page 121 *General settings*)


### Changing the sorting sequence of columns

- ◆ Click in the column header, drag and drop it at the required position.  
The columns are shifted according to the new position.

### Changing the column width

- 1 Move the mouse pointer onto the line in the header separating the columns.

- 2 Click and drag the column to the width you want.

### Adjusting to optimal width

- ◆ Double-click the column line in the title bar.

## 4.4.3 Searching for studies/patients

- ✓ Patient mode

### Updating the Pre-registered Patients list

If your system is connected to an HIS/RIS patients coming in are displayed in the **Pre-registered Patients list (worklist)**.

You can configure the update interval for the **Pre-registered Patients list**.  
DICOM properties

You can update the **Pre-registered Patients list** manually:

- ✓ **Pre-registered Patients list**
- ◆ Click this icon.

## Searching for studies/patients in a study list

- 1 Click the tab of the required study list.

- 2 Select the time period from the filter.
- 3 Enter any known **ID**:, **Name**:, or **AccNo**: (accession number) of the required patient if necessary.
- 4 Click **Search**.
- 5 To clear the search criteria: Click **Clear** and then **Search**.

## Searching for certain studies/patients from HIS/RIS

If the patient is already waiting you can ask the HIS/RIS for special entries only.

- ✓ The system is connected to an HIS/RIS.
- ✓ You have the permission to access the function.

- 1 Click this icon.

- 2 Enter any known data, if required.
  - ID
  - Name
  - Accession No.
  - Requested Procedure ID
  - Physician
  - SPS Description (Scheduled Procedure Step Description)
  - AE-Title
- 3 Click the check boxes of **Date** and **Time**, and narrow down the time period if necessary.
- 4 Click **Search**.

#### 4.4.4 Selecting a study/patient

For each study, one entry is assigned to in the study lists.

For each patient, there may exist one or more study entries.

If you want to use the data of a study for an examination, for postprocessing or documentation, or for a manual composing:

- 1 Click the correct entry in a study list.

The entry appears inverted.

- 2 Click this icon.

The images of the selected study are displayed.

Double click the required image (**Stitch** mode for an Ortho composed image and **View** mode for a normal image).

Click the other mode tab as required, and the system switches to the respective mode.

– or –

To switch to the default mode (**Examination for Pre-registered Patients, View for Examined Patients**):

Double-click the correct line in the study list.

You can now acquire new images for this study or stitch/postprocess/view/print images that already exist.

### 4.4.5 Patient registration

If your system is not connected to an HIS/RIS or in an emergency, you must manually register the patient to be examined.

Registration is necessary, since acquisition is not possible without it.

- ✓ Patient mode
- ✓ Pre-registered Patients list

### Emergency registration

In case of an emergency, you can perform a "fast registration" and enter the personal data of the patient later.

- 1 If you do not have a user account for this system: Click **Emergency** on the login screen.

– or –

If you are already logged in: Click this icon.

A new study for the emergency patient is registered.

The name of the emergency patient is set to: **ER\_**current date and time.

The new study is passed to the Examination mode.

- 2 Select the organ programs to be performed.

(→ Page 32 *Entering procedure information*)

- 3 Click **OK**.

Later on you can enter the missing personal data of the patient.

(→ Page 34 *Correcting study data*)

### Manual registration

For patients that are to be registered manually (not for HIS/RIS patients), you can enter the data as follows:

To create a new study/patient entry:

- 1 Click this icon.

– or –

To add a study for a patient already registered:

Select an item in the patient list and click this icon.

(Example configuration)

## 2 Enter all the relevant data.

Mandatory entry fields are indicated by asterisks \*.

You can configure the displayed items and the mandatory entry fields.

(→ Page 127 *Patient registration configuration*)

## 3 Select the procedures/organ programs to be performed.

(→ Page 32 *Entering procedure information*)

## 4 To register the patient and examine them later: Click **Pre-register**.

## 5 To clear the existing patient data: Click **Clear**.

## 6 To examine the patient immediately: Click **Examination**.


## 7 To return to Patient mode: Click the icon.

# Entering procedure information

Before starting the examination, you must select one or more procedures/organ programs.

(→ Page 163 *RIS mapping*)

### 1 Select the category.

- **Adult:** Procedures for adult
- **Paediatric:** Procedures for children

The age for pediatric patient can be set in the **Patient Registration Configuration**.

(→ Page 127 *Patient registration configuration*)

- **Often Used:** Favorite procedures


The favorite procedures can be set in the **Procedure data** in the OGP Editor.

(→ Page 144 *Procedure data*)

- **Stitch:** Procedures for Ortho

Before selecting a procedure for Ortho, you must select the patient type: **Adult** or **Paediatric**.

### 2 Select the body region by clicking in the figure.

### 3 Enter a search criteria for an organ program, if necessary.

- 4 Change the sort order for the procedures by clicking an icon, if necessary.

Custom sort order

You can define the customer sort order in **Examination Properties**.

(→ Page 124 *Examination properties*)

- or -

If you are authorized by your administrator, you can change the procedure order manually via drag and drop. Click the save button to adopt the custom sort order.

Ascending sort order

Descending sort order

Sort alphabetically

- 5 Double click a procedure.

The organ programs linked to the procedure is displayed in the OGP list.

- or -

Click the arrow down.

- 6 Repeat the steps for further organ programs.

To delete an organ program from the procedure:

- 7 Select the organ program and click this icon.

To delete all the organ programs from the procedure:

- 8 Click this icon.


#### 4.4.6 Correcting study data

If a patient name is misspelled, for example, you can correct the personal data or enter additional data of a patient already registered in the study list.

- 1 Click the corresponding entry in the study list.

- 2 Click this icon.

(Example configuration)

- 3 Select the study to correct, if necessary.
- 4 Correct or complete the personal data of the patient.

Patients provided by HIS/RIS can only be modified once their status has been switched to a local patient.

Once an HIS/RIS patient is switched to the **Examination** mode, they become a local patient. That means, the patient will further on not be regarded by HIS/RIS updates anymore.

- 5 Click **Accept**.

#### 4.4.7 Deleting studies

If pre-registered studies in the **Pre-registered Patients** list or (archived) studies in the **Examined Patients** list are not needed anymore, you can delete them.

You can configure automatic deletion of studies.

Disk management

You can also delete studies of a date range.

(→ Page 158 *Clearing study and image data*)

- 1 Click the required entry in a study list.

– or –

Press the **Shift** key or click the check box to select multiple entries.

The selected entries appear inverted.

- 2 Click this icon.

If there are studies with unarchived images among those selected, a further confirmation dialog appears.

- 3 If you want to delete the unarchived studies as well: Click Yes.

A confirmation dialog appears.

- 4 Check/change your selection.

- 5 To delete the selected entries: Click Yes.

When the **Recycle Bin** has been enabled, the deleted studies are moved to the **Recycle Bin**.

(→ Page 116 *Managing studies in the Recycle Bin*)

When the **Recycle Bin** has not been enabled, the studies are finally deleted.

#### 4.4.8 Protecting/unprotecting studies

You can protect studies from unintentional deletion.

- 1 Click the correct entry in a study list.

– or –

Press the **Shift** key or click the check boxes of the desired studies to select multiple entries.

The selected entries appear inverted.

To protect the studies:

- 2 Click this icon.

Protected studies are marked by a closed lock symbol.

To unprotect the studies:

- 3 Click this icon.


Unprotected studies are marked by an open lock symbol.

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
- 1 Connect the Core XL to the all-in-one PC with the cable.

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## Detaching the detector

- ✓ The detector has been registered before to this system by the Siemens Healthineers Service.
  - ✓ **Detector Overview** dialog is displayed.
- 1 Click **Disconnect** at the detector overview window.  
A message appears asking if you want to detach the detector.
  - 2 Click **OK** to start the detachment procedure.

A message appears stating that the detachment was successful.

The detached detector is removed from the detector overview window.

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A wireless detector can only be connected to your system if it has been registered for your system by its serial number and if it has been calibrated for use with your system.

Changing the detector is only possible in **Patient** mode.

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## Upgrade detector firmware

For optimum image quality, check and upgrade the detector firmware.

### ✓ FD Manager

✓ A new detector firmware version is available.

✓ Detector battery level > 50%.

#### 1 Click **Update**.

The upgrade is started.

A message appears stating that the firmware update is in progress.

Wait until the installation is done (approx 10 minutes).

A message appears stating that the firmware update is finished.


#### 2 Click **OK**.

## 4.4.10 Quality control study

There is a special study for daily quality control.

## Opening the quality control study

- ✓ You are logged in with a user account which has the **Daily Control** permission.
- ✓ **Patient mode**
- ✓ **Pre-registered Patients list**
- ◆ Click this icon.

You can now perform your daily quality control, acquire images, postprocess them, view the exam protocol, etc.

## Acquiring images for quality control

- 1 Click **Exam**.
  - 2 Select an organ program.
  - 3 Acquire quality control images.
- Acquisition data will be added in this list.

## Checking quality control images

- ◆ Double-click an entry in the list.

## Printing the quality control study log on a paper printer

- 1 Select the **Start Date** and **End Date**.
- 2 Click **Print**.

## Exporting the quality control study log

- 1 Plug-in the USB memory, if necessary.
- Wait a moment.

- 2 Select the **Start Date** and **End Date**.
- 3 Click **Export**.

The log file is saved in .csv format

### Deleting the quality control study log

- 1 Select the required entry.
- 2 Click **Clear**.

## 4.5 Image acquisition

- (→ Page 41 *General information*)
- (→ Page 45 *Checking/changing organ programs*)
- (→ Page 47 *Acquisition guidance*)
- (→ Page 48 *Checking/changing acquisition parameters*)
- (→ Page 54 *Organ program and acquisition parameters on the TUI*)
- (→ Page 55 *Releasing an exposure*)
- (→ Page 56 *Exposure index*)
- (→ Page 58 *Accepting or rejecting images*)
- (→ Page 60 *Repeating an exposure*)
- (→ Page 60 *Image recovery*)
- (→ Page 61 *Closing a study*)
- (→ Page 39 *Quality control study*)

Please also refer to the Operator Manual of the MULTIX Impact C.

### 4.5.1 General information

You can release acquisition in **Examination** mode.



If the imaging system is not ready, you cannot perform acquisitions.

Pay attention to error messages before releasing an acquisition.

(→ Page 19 *Status information*)

Patient data and system data are shown top left in **Examination** mode.

**The grid to be used for acquisition is defined in the organ program.**

**Grid type:**

- **Universal:** universal grid,  $f_0=140$  cm
  - **115 cm:** grid,  $f_0=115$  cm
  - **180 cm:** grid,  $f_0=180$  cm
- or -

### Grid type defined in the organ program

- **Current Grid: none:** No grid is inserted.
- **Current Grid: unknown:** For free exposures, a clip-on grid may be used. The status of the clip-on grid cannot be detected by the system.
- **Current Grid: removed:** Grid is removed during grid test. Table or wall mode is active, grid runs out of its position during grid test.
- **Desired Grid: XXX:** The grid type defined in the organ program is indicated.
- **Desired Grid: none:** No grid is defined in the organ program.
- **Desired Grid: unknown:** An organ program for free exposure has been selected. The grid is not defined in the free exposure organ programs.

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Please note that for perfect image quality, the grid focus ( $f_0$ ) must be within the admissible source image distance (SID) range!

Check the correct grid focus ( $f_0$ ) on the grid label before making acquisitions.

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#### Poor image quality leading to repeated X-ray exposure

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### Source to image distance (SID)

**SID:** The current SID is indicated.

For free exposures onto a detector or cassette, measure the SID with the tape measure in the multileaf collimator and adjust the SID value on the TUI.

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### Collimation

**Collimator Size:** The current collimation size from the automatic collimator is indicated.

For automatic collimator, the collimation size defined in the organ program is applied and indicated. You can change the collimation size by adjusting the two knobs on the collimator.

#### Collimator Alignment:

- **Center:** Center alignment is reached.
- **Center not reached:** Center alignment is not reached.
- **Top:** Top alignment is reached.
- **Top not reached:** Top alignment is not reached, although the function button on the TUI is pressed.

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## 4.5.2 Information on acquisition using a detector

(→ Page 21 *Detector status*)

### Detector calibration needed

If the flat detector needs calibration, a message is displayed and a calibration symbol is displayed in the status area.

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### Detector temperature

To obtain optimal image quality, the detector temperature must be within a certain range.

example: out of detector calibration  
temperature

If the temperature of the flat detector is out of calibration temperature range, a thermometer symbol is displayed in the status area.

(→ Page 21 *Detector status*)

Do not make acquisitions and wait until the detector temperature is within the required range for optimum image quality.

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**Poor detail resolution with all types of radiation!**

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### Detector battery charge

A battery symbol in the status area displays the charging condition of the selected detector.

### Detector wireless connection

A WIFI symbol in the status area displays the wireless signal strength of the selected detector.

### Preview image

The original detector images must be processed by some algorithms. To speed-up image display, a preview image is initially displayed.

After that, the fully processed image soon appears.

### 4.5.3 Checking/changing organ programs

During patient registration, organ programs are normally chosen for the planned procedure.

Examinations are based on organ programs. Each organ program can be assigned up to four workplaces: patient table, wall stand, free, and cassette.

The setup of the acquisition mode and workplace is done with the selection of an organ program.

When entering the **Examination** mode, the first unperformed organ program in the organ program schedule list for the active patient is automatically selected.

The selected organ programs are displayed in the **OGP List** or **Ortho OGP list**.

### Selecting an organ program for acquisition


- ✓ At least one organ program has been chosen during patient registration.
- ◆ Click the organ program in the **OGP List** to use.

The organ program is marked by a frame.

When adding a new procedure in the **Examination** mode, the new organ program will be saved in the study of the OGP which is currently marked by a frame in the **OGP List**.

(→ Page 46 *Adding organ programs*)

Example

— OR —

Double-click the **Ortho** organ program in the **Ortho OGP list** to switch the system to the **Ortho** examination mode.

(→ Page 65 *Organ program for Ortho*)

You can configure the sort order for organ programs.

(→ Page 124 *Examination properties*)

### Adding organ programs

- ◆ Click this icon.

(→ Page 32 *Entering procedure information*)

### Removing organ programs

- 1 Click the organ program to remove.  
The organ program is marked by a frame.
- 2 Click this icon.

### Repeating organ programs

- 1 Click the organ program to repeat.  
The organ program is marked by a frame.
- 2 Click this icon.

### 4.5.4 Acquisition guidance

Guidance information can be displayed for easy patient positioning.

- ✓ Examination mode
- ✓ The correct organ program has been selected.
- ◆ Click this icon.

An image indicating the patient positioning is displayed.

The position guide images are for reference only. The operator may position the patient according to his/her own knowledge. The correct positioning is the responsibility of the operator.

Not all organ programs contain position guide images.

On the **myExam 2D Camera** tab card, the live image from the camera in the examination room is shown (only when the camera is configured with the automatic collimator).

When positioning patient on the table, do not use the live image as position assistance if the SID is smaller than 115 cm.

On the **X-Ray Parameter** tab card, the default parameters for the selected organ program are shown.

#### 4.5.5 Checking/changing acquisition parameters


##### **CAUTION**

Wrong radiation parameters

**Risk of undesired X-ray exposure due to wrong parameters**

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✓ The correct organ program has been selected.

Example

## Workplaces and exposure techniques

The following combinations of workplace, acquisition mode, and exposure technique are possible:

| Icon | Workplace                                  | Acquisition mode         | Exposure technique |                                |         |         |
|------|--------------------------------------------|--------------------------|--------------------|--------------------------------|---------|---------|
|      |                                            |                          | 1-point            | 1-point with extended exposure | 2-point | 3-point |
|      | Tube with wireless FD in tray (table)      | DR RAD                   | x*                 | x*                             | x       | x       |
|      | Tube with wireless FD in tray (wall stand) | DR RAD                   | x*                 | x*                             | x       | x       |
|      | Tube with wireless FD (free positioned)    | DR RAD                   |                    |                                | x       | x       |
|      | Tube with cassette (free positioned)       | Cassette (free exposure) |                    |                                | x       | x       |

x\*: with automatic exposure control (AEC)

To change the workplace:

- ◆ Click the required icon.

## Patient size

The acquisition parameters can be adapted to the patient size.

Small/thin patient

Medium patient

Large/thick patient

To change the patient size:

- ◆ Click the required icon.

## Exposure technique

The following generator exposure techniques are supported:

- 1-point technique (with Automatic Exposure Control, AEC):  
An acquisition is taken with operator specified kV value.  
1-point technique with extended exposure (AEC with current reduction)
- 2-point technique:  
An acquisition is taken with operator specified kV and mAs values.
- 3-point technique:  
An acquisition is taken with operator specified kV, mAs and ms values.

Exposure time cannot exceed 450 ms for examinations with a wireless detector, or 5 s with a cassette.

The available exposure techniques depend on the system and the selected workplace.

The following table shows how you can change the exposure technique:

|                  |                   | Current technique                              |                                |                                                |                     |
|------------------|-------------------|------------------------------------------------|--------------------------------|------------------------------------------------|---------------------|
|                  |                   | 1-point                                        | 1-point with extended exposure | 2-point                                        | 3-point             |
| Target technique | 1-point           | -                                              | Deselect extended exposure.    | Click on measuring field icon or Dose field.   |                     |
|                  | 1-point ext. exp. | Select extended exposure.                      | -                              | Not possible - First select 1-point technique. |                     |
|                  | 2-point           | Click in mAs field.                            |                                | -                                              | Click in mAs field. |
|                  | 3-point           | Not possible - First select 2-point technique. |                                | Click in ms field.                             | -                   |

### Tube voltage kV

- ◆ Click + or -.

### Tube mAs

- ◆ Click + or -.

### Tube ms

- ◆ Click + or -.

### Dose (Sensitivity)

For acquisitions with automatic exposure control (1-point technique), you set the applied dose and therefore the sensitivity.

Adjusting the measuring fields to the region of interest greatly improves image quality.

When the dose is reduced (sensitivity is increased), the image quality might be sufficient for the type of examination and the region of interest.

The decision for a certain image quality and applied dose must be made by the physician or radiologist on the basis of experience, measured data, and the type of examination.

- ◆ Click + or -.

### Density (Exposure correction)

Exposure correction can be used to compensate the influence of e.g. metal implants.

The exposure correction display shows an initial value defined in the organ program.

This value can be manually modified in increments of  $\pm 0.5$  EP.

- ◆ Click + or -.

### Extended exposure (current reduction)

The current reduction extends the acquisition time for correct exposures in 1-point technique, for example, thorax of very thin patients or infants.

- ◆ Click this icon.

### Measuring field (dominant)

Depending on the organ to be examined, you must select one or more suitable measuring fields so that the automatic exposure control can match the dose to the organ of interest.

Left measuring field

Middle measuring field

Right measuring field

You can only select the left and/or right chamber if one or both of them are covered by the collimation.

- ◆ Click the required measuring field icons.

The measuring fields are indicated in the image.

During collimation, the used measuring field may be automatically changed to the middle measuring field.

Outlines of the active measuring fields are visible in the monitor image.

Pay attention to correct change-over.

When the camera is configured, the measuring fields are indicated at the upper right corner of the live image displayed in the **Position Assistance** tab card.

The indication of the measuring fields is for reference only and it does not correspond to the actual coverage of the measuring fields on the patient.

## Focus

You can either select the small or the large focus.

Small focus

Large focus

- ◆ Click the required focus icon.

## Tube load

You can select the tube load as either 80% or 100%.

When using tube utilization of 80%, exposures are made with smaller tube currents and longer exposure times. This way the X-ray tube is utilized in a safer way and will have a longer life.

100% tube utilization is nevertheless also admissible.

- ◆ Select the required value in the Tube Load: drop-down list.

## Copper filter

You can select one additional copper filter into the beam to harden the radiation, for example for pediatric examinations.

- ◆ For automatic collimator, select the required value in the Cu filter: drop-down list.

For more information, refer to the Operator Manual of the MULTIX Impact C.

## Top alignment

There are two collimation methods on the receptor stand:

- Center alignment: The X-ray radiation is collimated so that its central beam is centered to the detector.
- Top alignment: The X-ray radiation is collimated so that its top edge of the beam is aligned to the top edge of the detector.

Top alignment on the receptor stand is useful for acquisitions, e.g. for thorax examinations, where you first adjust the detector height to the patient size and then perform further collimation.

Top alignment is only possible for ACSS T (Automatic Collimation Size Sensing and Tracking) variants. For a manual collimator, top alignment is not possible.

To use top alignment:

- 1 Press this icon on the touch user interface.  
The light localizer is switched on.
- 2 Check the collimation.

### Resetting parameters to the organ program settings

If you want to reject parameters you have changed:

- ◆ Click this icon.

### Saving the organ program settings

The values of kV, mAs, ms, dose, density, measuring field, tube load, Cu filter, and focus can be changed and saved as default organ program.

If you are authorized by your administrator, you can save parameters you have changed as default organ program:

- 1 Click this icon.  
The following message appears: <"Are you sure you want to continue?">
- 2 Click OK.
- 3 If you have changed the exposure technique, the following message appears: <"You have changed the exposure mode of the current patient size, please make sure all three patient sizes use the same exposure mode.">
- 4 Click OK.


Always make sure that the same exposure technique and proper parameters have been set for all the three patient sizes.

#### 4.5.6 Organ program and acquisition parameters on the TUI

You can also select the organ program and check/change the acquisition parameters in the examination room on the touchscreen user interface (TUI).

For further functions refer to the Operator Manual of the MULTIX Impact C system.

#### Example

- ◆ Select the required values in the drop-down lists.

### 4.5.7 Releasing an exposure

Please also refer to the Operator Manual of the MULTIX Impact C.

- ✓ The correct flat detector is connected and selected, or a cassette has been placed (for free exposures only).
- ✓ **Examination** mode
- ✓ The patient is positioned correctly.
- ✓ The correct organ program has been selected.
- ✓ The correct workplace is selected.
- ✓ The correct acquisition parameters are set.
- ✓ The imaging system is ready for acquisition.

- 1 Press the acquisition push button.  
X-ray is released.

The image is processed and displayed in **Postprocessing** mode.

- 2 Check the image.  
(→ Page 76 *Image postprocessing*)
- 3 If the acquired image is okay: Click **Accept**.  
– or –  
If the acquired image is not okay: Click **Reject**.

If the image display is not satisfactory, you can correct it in postprocessing.  
(→ Page 76 *Image postprocessing*)

**CAUTION**

Image was displayed and then lost after exposure.

**Risk of repeated radiation**

- ◆ Call Siemens Healthineers service to restore the image.

You can configure automatic acceptance of each acquired image.

(→ Page 124 *Examination properties*)

After an acquisition has been performed, the next unused organ program is selected automatically.

(→ Page 7 *Image failures*)

#### 4.5.8 Exposure index

To avoid systematic errors, like under- or overexposure, the imaging system provides the Exposure Index (EXI) as an initial indication for the applied dose.

The exposure index is an important indicator for the exposure of the acquisitions.

The exposure index is only intended to be a coarse indicator of the dose applied to the receptor for a particular examination. The magnitude of this unitless value depends on many factors. It is therefore a relative value under certain acquisition conditions.

There are several definitions for the exposure index:

- The Clinical EXI analyzes the mean gray value of the diagnostic relevant part in the image via histogram, without the peaks for direct radiation and collimation.
- The Deviation Index detects the tolerance of the exposure in correlation of the Clinical EXI value. The Deviation Index tells something about under- and overexposure of the diagnostic area.

As a guidance, a target clinical EXI value can be defined in every organ program.

The Clinical EXI value is displayed on the bottom-left of the image.

- The target clinical EXI value (**Reference EXI:**) which is defined in the organ program (**Target EXI**).

(→ Page 153 *Miscellaneous parameters*)

- The Clinical EXI value of the current image (**Current EXI:**)
- The Deviation Index of the current image (**DI:**)

The exposure index values, and the Deviation Index are stored as private tags and as (new) public tags in the DICOM header of an image. For details, see the DICOM Conformance Statement.

The exposure index values are therefore included in image transfers.

The exposure index values are also included in the exam protocol.

(→ Page 102 *Displaying the exam protocol*)


Pay attention to the exposure index in all examinations.

## Clinical EXI

The Clinical EXI or EI is determined as the "value of Interest" of the image. The calculation considers the relevant area of the image and eliminates non relevant areas like directly radiated areas or non-collimated areas.

The Clinical EXI considers the image content.

The Clinical EXI is normalized to calibration conditions, i.e. it is widely independent from systematic influences, e.g. from the used flat detector. For this purpose, a normalization factor can be entered in service configuration.

## Linear relation

Within defined conditions, the Clinical EXI is directly proportional to the applied dose (linear relation), i.e. a double in the exposure index corresponds to twice the dose.

If the same organ program is selected under similar acquisition conditions and the same positioning of the patient has been maintained, the exposure index must be within a certain range.

Example:

Thorax series examination

With identical positioning, e.g. central p.a., and exposure parameters, the exposure index value must only fluctuate within certain limits.

Please note, that the reference value for the exposure index depends from the examination type and the organ examined. The exposure index values differ also with different projections/views.

## Dependencies

The exposure index depends on the following factors:

- selected organ program
- organ examined
- projection

- collimation
- position of the organ examined
- metal implants, pacemakers, inside the considered area
- position of the AEC dominant (measuring field)

## Large deviations

If the exposure index deviates greatly from the usual values, this might be due to incorrect positioning, incorrect setting, or in extreme cases, a defect in the system.

If the exposure index deviates greatly from the usual values during a standard examination and other causes for the deviation have been eliminated:

Perform the constancy test.

If this test reveals deviations from previous measurement results:

Call Siemens Healthineers Service to have the system inspected.

If no deviations are found, the observed deviation was probably due to the dependencies described above.

## Deviation Index

The Deviation Index or DI is a gauge for the deviation of the Clinical EXI of an image from the target value scaled in exposure points.

The Deviation Index is calculated after image acquisition by the following formula:


$$DI = 10 \cdot \log_{10} (EI / EI_T)$$

EI is the clinical exposure index of the actual exposure (**Current EXI**).

$EI_T$  is the target clinical exposure index (**Reference EXI**) for the examination, which is defined in the organ program (**Target EXI**).

(→ Page 153 *Miscellaneous parameters*)

### 4.5.9 Accepting or rejecting images

Example

Accepted image

Rejected image

You can configure automatic sending/archiving of images.

DICOM properties

The images will be sent automatically to the configured sending/archiving nodes.

Only accepted images will be sent.

✓ **Examination mode**

✓ One or more images have been acquired.


### Rejecting an image

1 Click the image to reject.

The image is marked by a frame.

2 Click this icon.

Select or enter a reject reason and click **OK**.

You can configure the predefined reject reasons.

(→ Page 122 *Image processing properties*)

### Accepting a rejected image

1 Click the image to accept.

The image is marked by a frame.

2 Click this icon.

### 4.5.10 Repeating an exposure

You can reselect a used (checked) organ program that has been previously used for acquisition.

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✓ An exposure has been rejected.

1 Click the image/organ program to repeat.

The image/organ program is marked by a frame.

2 Click this icon.

3 (→ Page 55 *Releasing an exposure*)


### 4.5.11 Image recovery

If an image transfer could not be completed from the flat detector, e.g. due to a connection problem, the acquired images will not be lost. The images from the flat detector will be transferred at a later time, when connection is up again.

In case of a battery failure the images are lost, an image recovery will not be possible anymore.

A window appears on the monitor.

1 Make sure that the detector is close to the access point.

If connection can be established, the images are transferred to the imaging system.

- 2 Click **OK** and wait until the window disappears.

A window appears if the image recovery fails after a timeout.

- 3 Click **Retry** to try again.

If you click **Cancel** to abort the image recovery, the image will be lost.

If image recovery is aborted:

- ◆ You may continue working by reconnecting the wireless detector.
- ◆ Reset the detector. Do so by pressing the power button (for MAX wi-D) until all the LED indicators are off or by taking out the battery.

In this case, the image will be lost.

After switching the detector on again or reinserting the battery, you can resume working with the detector. Make sure you have a good and continuous WiFi connection when working with the detector.

#### 4.5.12 Closing a study

When an examination is finished, the study is either closed automatically or you can close the study manually.

You can configure that the study close dialog pops up automatically when reentering the **Patient** mode with an open study.

(→ Page 176 *Configuring the MPPS node*)

#DROC DICOM properties#

### Automatic MPPS close

- ✓ Automatic MPPS close is configured.

- ◆ Click the **Patient** tab.

The study is closed and stored in the **Examined Patients** list if all organ programs of the study have been performed.

Otherwise, the **MPPS Dialog** is displayed.

(→ Page 63 *MPPS close*)

### Manual MPPS close

- ✓ Manual MPPS is configured.

- 1 Click the **Patient** tab.

The **MPPS Dialog** is displayed.

(→ Page 63 *MPPS close*)

- 2 Select the **Status**.

**In Progress:** The study is still in progress and stored in the **Pre-registered Patients** list.

**Completed:** The study is closed and stored in the **Examined Patients** list.

**Discontinued:** The study is aborted and stored in the **Examined Patients** list.

- 3 Click **OK**.

Studies set as **Completed** or **Discontinued** cannot be reopened for acquisition.

In the **Examined Patients** list, the studies set as **Completed** are marked with this

icon

### No MPPS configured

- ✓ MPPS node is disabled. (→ Page 177 *Configuring sending/archiving nodes*)

- ◆ Click the **Patient** tab.

The study is set as **Completed** and closed if all organ programs of the study have been performed.

Otherwise, the study is set as **In Progress**.

You can close a study in the **Pre-registered Patients** list by clicking this icon


- If MPPS is configured, the **MPPS Dialog** is displayed.
- If MPPS is not configured, the selected study is set as **Completed** and stored in **Examined Patients** list.

#### 4.5.13 MPPS close

If your institution uses MPPS (Modality Performed Procedure Step), an examination of a patient has to be confirmed when interrupted or discontinued.

In this case, the **MPPS Dialog** will be displayed.

- **Patient:** Patient data is displayed.
- **Dose:** Dose data is displayed.
- **Films:** Used film quantities.
- **Study:** Study data is displayed.
- **Performed Procedure:** Performed procedure step information is displayed.
- **Performed Protocol Codes:** Select the performed protocol code, if necessary.
- **Status:** Select the study status.
  - **In Progress:** The examination is still in progress and stored in the **Pre-registered Patients** list.
  - **Completed:** The study is closed and stored in the **Examined Patients** list.
  - **Discontinued:** The study is aborted and stored in the **Examined Patients** list. Select a reason for the cancellation.

A study marked as **Completed** or **Discontinued** cannot be continued.  
Auto sending or filming will not be performed for the study marked as **Discontinued**.

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(→ Page 65 Organ program for Ortho)

(→ Page 65 Preparing Ortho)

(→ Page 70 Releasing the Ortho exposures)

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
### Preparing Ortho

Pay attention to the messages displayed on the imaging system display and on the TUI. If the series acquisition is interrupted, the whole sequence must be started again from the beginning!

(→ Page 65 *Organ program for Ortho*)

- ✓ An Ortho organ program is selected.
- ✓ The collimator must be in the 0° lock-in-position.

- ✓ For Ortho acquisitions onto the wall stand, the detector unit must be in the 0° position.

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(→ Page 66 *Positioning the patient on the patient table*)

(→ Page 66 *Positioning the patient on the Ortho stand*)

#### Positioning the patient on the patient table

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#### **Risk of crushing between table and detector tray**

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#### Positioning the patient on the Ortho stand

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- 2 Position the patient standing on the Ortho stand in such a way that the longitudinal axis of the object runs parallel to the image receptor.

- 3 Position the calibration element (sphere, ruler) in the area relevant for diagnosis.

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- ✓ After selecting an organ program for Ortho, the TOD defined in the organ program is displayed.

- ✓ The X-ray tube unit is positioned for Ortho.

- 1 Measure the TOD using a ruler.

– or –

Measure the SOD (source-object distance) using the tape measure.

(TOD = SID - SOD)

- 2 Select the correct TOD on the TUI or imaging system (accuracy 0.5 cm).

#### Adjusting the ROI with start and end position in the examination room

- ✓ SmartOrtho is configured.

- ✓ (→ Page 65 *Preparing Ortho*)

#### Example of Ortho menu

| Icon | Short description                                                      |
|------|------------------------------------------------------------------------|
|      | Select TOD (Table-Object Distance)                                     |
|      | Acquisition start or end position (based on the acquisition direction) |
|      | Acquisition end or start position (based on the acquisition direction) |

- ✓ The Ortho menu is displayed on the TUI.

Use the laser line light localizer to set the image area.

- 1 Ask the patient to close their eyes and turn on the light localizer.
- 2 Align the laser line on the upper or lower edge of the region of interest by tilting the X-ray tube unit.

- 3 Press one of these two icons on the TUI.

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- 4 Align the laser line on the lower or upper edge of the region of interest by tilting the X-ray tube unit again.

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#### Adjusting the ROI with camera in the control room

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- ✓ Smart Virtual Ortho is configured.
- ✓ (→ Page 65 *Preparing Ortho*)
- ✓ The camera is configured for Ortho.
- ✓ The X-ray tube unit is positioned for Ortho.
- ✓ TOD has been set correctly.

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#### ROI auto detection

- ✓ Smart Virtual Ortho is configured.
- ✓ Auto Long-Leg/Full-Spine Collimation is configured.
- ✓ (→ Page 65 *Preparing Ortho*)
- ✓ The camera is configured for Ortho.
- ✓ Workplace wall stand is selected.
- ✓ The X-ray tube unit is positioned for Ortho.
- ✓ TOD has been set correctly.

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- You can manually adjust the ROI.
- Use the additional removable platform on the Ortho stand when performing Ortho acquisition on the wall stand.
- Adjust the upper or lower offset in the **Ortho Setting** function.

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#### Accepting or rejecting the acquired images

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## 4.6 Stitch mode


In **Stitch** mode, you can view and perform manual composing for the Ortho series.

The following functions are available:

- **View:** Functions for changing the image display and performing measurement
- **Composing:** Functions for manual composing

### 4.6.1 Viewing an Ortho image

The **View** tools provide various functions for changing the image display and performing measurements.

#### Example

- (→ Page 83 *Panning an image*)
- (→ Page 82 *Zooming an image*)
- (→ Page 96 *Measuring angles*)
- (→ Page 96 *Measuring distances*)
- (→ Page 97 *Deleting measurements*)
- (→ Page 81 *Cropping an image*)
- (→ Page 72 *Adding a laterality label*)
- (→ Page 73 *Saving changes*)

### Adding a laterality label

You can add a laterality label (R/L label) and position it in the image.

You can configure the laterality label language.

(→ Page 125 *Image display properties*)

#### To add a left or right label

1 Click this icon to enable the label function.

Other functions in the **Tool** and **Measurement** panels are disabled.

Click the icon again to disable the label function.

2 Click one of these icons.


3 Drag and drop the label to the required position in the image.

To change the label size

✓ A laterality label has been placed in the image.

1 Click this icon.

2 Click the + or - icons or enter the required text size.

3 Click the icon again.

The laterality label in the image is changed.

To delete the laterality label

✓ A laterality label has been placed in the image.

◆ Click this icon.

The laterality label is removed.

## Saving changes

After loading a previously acquired Ortho series, you can save the postprocessing you have carried out on the image.

✓ An Ortho series has been processed.

◆ To save the image: Click **Save**.

### 4.6.2 Composing

After the Ortho acquisition is finished, the single images can be stitched together and displayed as an overlapping single image via **Composing**.

Composing will be performed automatically and can be adapted manually afterwards.

- Automatic Composing: Images will be composed automatically. The algorithm for image composition (spine or leg) is automatically selected.
- Manual Composing: You can also display an image from the series to compose in **Stitch** mode and then start composition manually.

#### CAUTION

Motion between Ortho images

##### Incorrect diagnosis basis

- ◆ Check results of composing for motion artifacts due to mismatched Ortho images.
- ◆ Compare the positions of the collimator blades in each image.
- ◆ If the blades are correctly matched, but the anatomy shows a mismatch, then the patient moved between images.

If you are not satisfied with the auto composed image, you can perform the manual composing with the acquired Ortho series.

(→ Page 74 *Performing Composing manually*)

You can postprocess and document composed images in the normal manner. However, not all postprocessing functions are available.

## Performing Composing manually

You can manually compose an Ortho series which has just been acquired or which is stored in the local data.


1 Select the Ortho series to compose in the **Examined Patients** list.

2 Double-click the patient.

The images of the Ortho series and the composed image are displayed in **View** mode.

3 Click the **Stitch** task card.

4 Switch to the **Composing** menu in the control area.

### Example

5 Select the images to compose.

6 Click **Start**.

The images are composed.

The composing result of the first two images is displayed.

7 Click **Next** to check the other part of the composed image until the last two images have been stitched..

8 Click **Save** to save the result.


– or –

Click **Back** to adjust the image transitions manually.

(→ Page 75 *Adjusting the image transitions manually*)


## Changing the display of a composed image

- ◆ You may change the display of the image transitions by changing the **Transition Mode** via clicking the **Blend** or **Transparent** button.

**Blend mode:** The single images are blended at the transitions between each pair of images.

(1) Interpolated areas

**Transparent mode:** The single images are placed side by side without blended transitions.

When **Blend** mode is used, there are transition areas between the single images that were interpolated from two images.

Although the **Blend** mode image might appear more pleasing to the eye, the diagnosis of any anatomy should only be based on the original images.

## Adjusting the image transitions manually

In rare cases, where image transitions have not been optimally detected, you can manually adjust the position of the images relative to one another (shifting pixels).

Proceed by shifting the image part below the transition edge relative to the image part above the transition edge until the best alignment is achieved in the overlapping area.

As the user, you are responsible for correct manual adjustment of the image transitions and the conclusions for diagnosis drawn from it.


If any manual adjustment was performed, the text “manual adjustments performed” is displayed in the image.

1 Adjust the view to the image transition to be changed.

(→ Page 75 *Changing the display of a composed image*)

2 Change the required number of pixels to be shifted per step from the **Steps** list.

- 3 Click the arrow once or repeatedly in order to shift the image below in the corresponding direction.

The shift is performed after each click or continuously until the arrow is released.

The topmost image cannot be shifted.

The shift should be performed in **Transparent** mode because the pixel shift is easier to see in this mode.

(→ Page 75 *Changing the display of a composed image*)

- 4 Click **Next** to display and adjust the other images in the composed image.

- 5 Click **Save** to save the result.

## 4.7 Image postprocessing

In **Postprocessing** mode, you can process and annotate your images with the image processing functions.

The following functions are available for postprocessing:

- **Process:** Functions for changing the image display and annotating  
(→ Page 76 *Processing images*)
- **Flavor Setting** and **Rad:** For dedicated organ groups, the imaging system provides several sets of pre-defined image impressions (flavors) based on postprocessing parameters to address different preferences on image quality and presentation.  
(→ Page 88 *Setting an image flavor*)  
(→ Page 92 *Using the Rad tool*)

### CAUTION

Wrong image processing parameters

**Risk of error in diagnosis/artificial artifacts**

- ◆ Apply processing parameters carefully, e.g. windowing.
- ◆ Always check image quality visually before performing a diagnosis.

### 4.7.1 Processing images

The **Process** tools provide various functions for changing the image display and annotating.

- (→ Page 78 Adding a laterality label)
- (→ Page 78 Adding texts)
- (→ Page 79 Deleting text items)
- (→ Page 80 Windowing an image)
- (→ Page 81 Inverting the gray scale of an image)
- (→ Page 81 Cropping an image)
- (→ Page 82 Zooming an image)
- (→ Page 83 Panning an image)
- (→ Page 84 Flipping an image)
- (→ Page 84 Rotating an image)
- (→ Page 85 Setting the position and laterality information)
- (→ Page 87 Undoing changes)
- (→ Page 87 Saving changes)

### 4.7.2 Adding a laterality label

You can add a laterality label (R/L label) and position it in the image.

You can configure the laterality label language.

(→ Page 125 *Image display properties*)

To add/change a left or right label:

- 1 Click one of these icons.
- 2 Drag and drop the label to the required position in the image.

(→ Page 79 *Deleting text items*)

### 4.7.3 Adding texts

You can add predefined texts or free texts and position them in the image.


You can configure the predefined texts.

(→ Page 122 *Image processing properties*)

- 1 Click this icon.

- 2 Select or enter a text item and click OK.
- 3 Drag and drop the text to the required position in the image.


(→ Page 79 *Deleting text items*)

### 4.7.4 Changing the text size

- ✓ A laterality label or texts have been placed in the image.

- 1 Click this icon.

- 2 Click the + or - icons or enter the required text size.
- 3 Click the icon again.

The text size of all texts and the laterality label in the image is changed.

#### 4.7.5 Deleting text items

- ✓ A laterality label or texts have been placed in the image.

- 1 Click this icon.
- 2 Click the text item or laterality label in the image.  
A confirmation dialog is displayed.
- 3 Click Yes.


You can also right click the text item or laterality label to delete it.

#### 4.7.6 Windowing principle

The use of contrast and brightness windowing allows you to adapt the image to display the region of interest with optimal contrast conditions.

The individual pixels of an X-ray image are stored as integer values. These values represent a range of image intensity divided into increments much finer than the human eye can detect. The monitor can display a range of intensity values. The purpose of the windowing function is to determine a range of values to be displayed and the way they are adapted to be displayed on the monitor.

The displayed window area is defined by a pair of numbers. The values used for non-subtracted images are the Window Center = WC and the Window Width = WW.

Principle of windowing (display contrast agent/bones black)

- (1) Stored value range
- (2) Screen value range
- (3) WW: Window Width
- (4) WC: Window Center
- (5) Lowest absorption
- (6) Highest absorption
- (7) Dark value
- (8) Bright value

#### 4.7.7 Windowing an image

The current window values are shown in the bottom right-hand corner of the image.

|    |                                                                     |
|----|---------------------------------------------------------------------|
| WC | Brightness/Window Center<br>Value range 0 to 4095 for 12-bit images |
| WW | Contrast/Window width<br>Value range 1 to 4095 for 12-bit images    |

(→ Page 79 Windowing principle)

Extreme window values may reduce the image quality to a point where diagnosis is no longer possible.

Always check image quality visually before performing a diagnosis.

- ✓ An image has been acquired or loaded for postprocessing.
- ✓ Postprocessing mode, Process tools

- ◆ Press the middle mouse button (mouse wheel) on the image and move the mouse.
  - Upward: Decrease window center/brightness values.
  - Downward: Increase window center/brightness values.
  - To the right: Increase window width, decrease contrast.
  - To the left: Decrease window width, increase contrast.

– or –

Click an icon.

Increase window center/brightness values.

Decrease window center/brightness values.

Decrease window width, increase contrast.

Increase window width, decrease contrast.


The window values and current EXI are updated according to the current cropping mask.

#### 4.7.8 Inverting the gray scale of an image

- ◆ Click this icon.

The gray scale of the image is inverted (black/white inversion: bones black to bone white or vice versa).

#### 4.7.9 Cropping an image

Images are normally cropped automatically after acquisition.

Depending on the acquisition system, different cropping methods are applied:

- Auto cropping by using the collimation detection algorithm
- Auto cropping by using the collimator values
- Manual cropping

You perform cropping when manual cropping is required or if you want to change the automatic cropping.

#### CAUTION

Image cropping

**Image data is lost.**

- ◆ Check that the parts of the image relevant for diagnosis are still visible before storing a cropped image.

To remove the cropping mask: Click this icon.

To crop a rectangular region: Click this icon and drag a region in the image to redefine the shutter area.

To crop a rectangular region for a film size: Click this icon and select a film size from the list to redefine the shutter area.

To add a cropping mask: Click this icon, click the corners of the polygon region in the image, and finish by a double-click.

To fit the shutter to window: Click this icon.

The image is adjusted and centered to the cropping region and fit to the image display area.

Shutter is the display area of the image and displayed with a blue frame. Cropping mask is the processing area and defined by the result of cropping with a white frame. After acquisition and without any manual cropping, the cropping mask is the size of actual collimation and the shutter is the size of image presented.

Image area outside the shutter can not be sent, exported or printed.

You can configure defaults for the shutter and cropping mask.

(→ Page 122 *Image processing properties*)

### 4.7.10 Zooming an image

You can change the displayed size of the image by various functions.


### Using a magnifying glass

- ◆ Click this icon.

Click the detail in the image that you want to enlarge.

The image around the position of the mouse pointer appears enlarged by a factor of 2.

### Zooming the image continuously

- ◆ Click this icon.

Drag the mouse cursor in the image up or down.

– or –

Rotate the mouse wheel.

The image is resized simultaneously.

### Adjusting the image to the image display area

- ◆ Click this icon.

The image is resized so that the (cropped) image fits to the image display area).

### Adjusting the image to 1:1 pixel display

- ◆ Click this icon.

The image is resized so that each monitor pixel corresponds to one image pixel.

#### 4.7.11 Panning an image

- ✓ The image has been cropped or zoomed.

- 1 Click this icon.
- 2 Press the left mouse button and move the mouse cursor in the image.  
The image is panned simultaneously in both axes (behind the shutter).

#### 4.7.12 Flipping an image

The default image flip is defined in the organ program.

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##### Risk of incorrect diagnosis or treatment because of left/right or up/down confusion

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- 1 Click one of these or both icons.

The image is flipped according to flip axes that correspond to the vertical/horizontal center axes of the monitor (left and right / top and bottom are swapped round).

The flip indicator is displayed in the image.

The default position of the flip indicator depends on the settings in the organ program.

(→ Page 153 *Miscellaneous parameters*)

If a R/L laterality label is placed in the image, the label will be flipped accordingly.

To reposition the flip indicator:

- 2 Drag and drop the flip indicator to the required position in the image.

#### 4.7.13 Rotating an image

You can rotate an image by  $\pm 90$  degrees or by an arbitrary angle.

### Rotating the image by $\pm 90$ degrees

- ◆ Click one of these icons.

### Rotating the image continuously

- 1 Click this icon.

- 2 Enter the rotation angle or click +/-.
- 3 Click the icon again.


## 4.7.14 Setting the position and laterality information

The view position and image laterality have been defined in the organ program.

If the information is missing or not correct, you can change it.

- 1 Tap the icon in the **Static Tools** area.

- 2 Select the correct **Patient Position**, **View Position**, and **Image Laterality**, if required.
- 3 Confirm with **OK**.
- 4 Click **Save**.

## Patient position, view position, and image laterality

The Patient Position can be set as follows:

| Patient position             | Abbreviation | Symbol                                                                                |
|------------------------------|--------------|---------------------------------------------------------------------------------------|
| Head First - Supine          | HFS          |    |
| Head First - Prone           | HFP          |    |
| Head First - Decubitus Left  | HFDL         |    |
| Head First - Decubitus Right | HFDR         |    |
| Feet First - Supine          | FFS          |    |
| Feet First - Prone           | FFP          |    |
| Feet First - Decubitus Left  | FFDL         |   |
| Feet First - Decubitus Right | FFDR         |  |
| not defined                  |              |  |

The patient position entered affects the orientation of the image on the screen, i.e., the images are displayed as if the patient was standing in front of the viewer and facing the viewer.

The View Position can be set as follows (depending on examination):

- **AP:** Anterior-Posterior view
- **PA:** Posterior-Anterior view
- **lateral:** Lateral view
- **oblique:** Oblique view
- **not defined:** No view is set

Laterality information (Right R or Left L) is used to label paired organs, e.g. hands.

The Image Laterality can be set as follows:

- **right/R:** For right organ
- **left/L:** For left organ
- **both (R+L):** Paired organs, both shown
- **unpaired:** Unpaired organ
- **not defined:** No laterality is set

Laterality (right/left assignment) applies to the image content (right/left organ) and not to the right/left side of the patient.

### 4.7.15 Undoing changes

You can undo any image postprocessing you have made.

- ✓ Postprocessing mode, Process tools

#### Undoing the last change

You can undo a maximum of previous 10 steps.

- ◆ Click this icon.

#### Restoring the state of the image

- ◆ Click this icon.

The state of the image initially after loading is restored.

(→ Page 87 *Saving changes*)

### 4.7.16 Saving changes

- ✓ Postprocessing mode
- ✓ An image has been acquired or processed.

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
- ◆ If the image is okay: Click **Accept**.

– or –

If the image is not okay: Click **Reject**.

Select or enter a reject reason and click **OK**.

You can configure the predefined reject reasons.

(→ Page 122 *Image processing properties*)

## Saving an image

After loading a previously acquired image, you can save the postprocessing you have made to the image or save a copy of the image.

- ◆ To save the image: Click **Save**.

– or –

To save a copy of the image: Click **Save As**.

(→ Page 87 *Undoing changes*)

### 4.7.17 Setting an image flavor

The system provides four sets of pre-defined image impressions (flavors) for each organ program: standard, smooth, high contrast, and finely detailed.


After each acquisition, the image displaying on the screen is the one processed by default flavor parameters.

You can select one of the four flavors.

You can also define one customized flavor for each organ program.


- 1 Click **Flavor Setting**.

A confirmation dialog is displayed if you have changed the window value without saving in **Process** tab or current image is using a flavor other than the default flavor.

**2 Click Yes.**

The image processed with default flavor parameters is displayed.

**3 To change the flavor of the image, click the flavor icon that you want for the image.**

The image impression is changed to the chosen flavor.

**4 You may adjust the windowing of the image by dragging the **WCShift** (window center) or **WWShift** (window width) slider or clicking an arrow.**

(→ Page 87 *Saving changes*)

### Default flavor

To change default flavor for an organ program:


- 1 Click the flavor icon that you want for the organ program.
- 2 Click **Save as default**.

#### 4.7.18 Customized flavor

You may define one customized flavor for each organ program.

- 1 Click **Rad**.

A confirmation dialog is displayed if you have changed the window value without saving in **Process** tab or current image is using a flavor other than the default flavor.

- 2 Click **Yes**.

**3 Select Customized from the Set Property list.**

Adjust the parameters until the image reaches your requirements.

**4 Drag a slider or click an arrow.**

|                 |                                                                                                 |
|-----------------|-------------------------------------------------------------------------------------------------|
| <b>Contrast</b> | Increasing the value amplifies the image therefore the difference of grey levels is enhanced.   |
| <b>Detail</b>   | Increasing the value makes fine structures more visible.                                        |
| <b>Latitude</b> | Increasing the value reduces dynamic range and makes more grey levels visible at the same time. |
| <b>Noise</b>    | Increasing the value suppresses more noise and makes the image smoother.                        |
| <b>WC Shift</b> | Increasing the value increases the brightness of the image.                                     |
| <b>WW Shift</b> | Increasing the value decreases the contrast of the image.                                       |

**5 Click Save Params.**

The parameter values are saved into the customized flavor setting.

### Saving a customized flavor as default

If you are authorized by your administrator, you can save the customized flavor as the default flavor for the organ program:

- ◆ Click **Flavor Setting** and then click **Save as Default**.

#### 4.7.19 Using the Rad tool

You can adjust the image impression by using the **Rad** tool and save the adjusted results for the current image.

#### CAUTION

Reduced image quality due to filtering

**Risk of artifacts, such as metal implants appearing to be loose**


- ◆ Choose the latitude, noise, detail, and contrast carefully.

- ✓ An image has been acquired or loaded for postprocessing.

- ✓ **Postprocessing** mode


- 1 Click **Rad**.

A confirmation dialog is displayed if you have changed the window value without saving in **Process** tab or current image is using a flavor other than the default flavor.

- 2 Click **Yes**.

3 Select a Set Property.

4 Select **SIGMOID** in the LUT Curve Type.

**LINEAR** is used for quality check.

Adjust the parameters until the image reaches your requirements.

5 Drag a slider or click an arrow.

(→ Page 87 *Saving changes*)

Only parameter value changes for the **Customized** flavor can be saved into OGP database.

All parameter value changes for flavor **Standard**, **Smooth**, **High Contrast**, and **Finely Detailed** cannot be saved. The effect of adjustments only applies to the current image.


## Reset parameters

To reset the parameter value to the initial status:

- ◆ Click Reset Params.

## 4.8 Image evaluation

In **View** mode, you can view your images, perform measurements, select images for sending, and send them.

- (→ Page 95 Viewing images)
- (→ Page 82 Zooming an image)
- (→ Page 83 Panning an image)
- (→ Page 96 Measuring angles)
- (→ Page 96 Measuring distances)
- (→ Page 97 Deleting measurements)
- (→ Page 97 Calibrating images)
- (→ Page 100 Drawing arrows)
- (→ Page 101 Storing an image)

- (→ Page 101 *Marking images*)
- (→ Page 102 *Displaying the exam protocol*)
- (→ Page 102 *Sending/archiving images*)


### 4.8.1 Viewing images

- ✓ A study has been loaded for viewing the images.
- ✓ View mode

#### Displaying DICOM information

- ◆ Click this icon.

#### Changing the image display layout

- ◆ Click one of these icons.

#### Paging images

To display the previous/next image:

- ◆ Click one of these icons.

– or –

Use the scroll bar.

– or –

Rotate the mouse wheel.

### 4.8.2 Measuring angles

- 1 Click this icon.
- 2 Click the vertex point of the angle in the image.
- 3 Click the end point of the first line in the image.

Right click to restart the measurement.

- 4 Click the end point of the second line in the image.

(→ Page 97 *Deleting measurements*)

### 4.8.3 Measuring distances

The distance measured and showed on the display does not match the real geometric size of the object as no calibration has taken place.

You can recalibrate the image and/or measure distances.

All distance entries are automatically recalculated when you recalibrate the distance measurement.

#### CAUTION

Measurements made on an image that is not true-size

#### Incorrect diagnosis possible

- ◆ Do not assume that images are displayed in true-size on the monitor or film.
- ◆ Use the distance function for critical measurements.

#### CAUTION

Misinterpretation of the image size on the film

#### Wrong measurements used for further treatment (e.g. prosthesis selection)

- ◆ Check the scale shown on the film.
- ◆ Calculate any distances measured on the film to correspond with this scale.

Distortions are possible with projection images.

To draw a distance line:

- 1 Click this icon.
- 2 Click the first point of the distance line in the image.
- 3 Click the second point of the distance line in the image.

(→ Page 97 *Deleting measurements*)

#### 4.8.4 Deleting measurements

- ✓ One or more measurements are drawn in the image.

##### Deleting a single measurement

- ◆ Click a measurement line in the image and then click this icon.

##### Deleting all the measurements

- ◆ Click this icon.

#### 4.8.5 Calibrating images

You can (re-)calibrate an image via geometry calibration or by drawing a known distance in the image.

Measurements created before image calibration will be cleared.

(→ Page 98 *Geometrical effects*)

#### CAUTION

Incorrect calibration

Incorrect diagnosis due to inaccurate distances

- ◆ The examiner is responsible for applying the correct calibration method.
- ◆ Make sure that the calibration object and the measured object are in the same plane.
- ◆ Perform calibration carefully.

- 1 Click the ... icon and select **Measurement Calibration**.

#### Measurement Calibration

Geometry

Fiducial

#### Tips:

If you ..., please select Geometry to calibrate.

If you ..., please select Fiducial to calibrate.

Cancel


- 2 Click **Geometry** or **Fiducial**.

(→ Page 98 *Calibrating by geometry*)

(→ Page 100 *Calibrating by a known distance*)

## Geometrical effects

Because of geometric relationships in images obtained from point-source x-ray exposures, object with different distances to the image receptor will be magnified differently. This parallax phenomenon is demonstrated in the diagram below.

Although the objects "A" and "B" are equal in size, they are magnified differently.

The parallax phenomena has more effect close to the image borders and will be smaller, the greater the distance between focus and object.

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
## Calibrating by geometry

Geometry calibration is performed by using the object height.

Geometry for distance calibrations (left with table, right without table)

- (1) Source plane (focus)
- (2) Collimator blades
- (3) Object plane (organ of interest)
- (4) Tabletop surface plane
- (5) Detector housing
- (6) Detector surface plane
- (7) X-ray tube
- (8) Collimator
- (9) Tabletop
- (10) Flat detector

- SOD: Source-Object Distance

- SID: Source-Image Distance

- STD: Source-Table Distance


- TOD: Table-Object Distance

- OS: Object size

- IS: Image size

- DS: Detector size

✓ (+ Page 97 Calibrating images)

- 1 Check/enter the Source-Image Distance **SID** (cm).
- 2 Check/enter the Table-Object Distance **TOD** (cm).
- 3 Click **Calibration**.

### Calibrating by a known distance


You can calibrate an image by drawing a distance line which measures the dimensions of a calibration object, e.g. a ruler.


The calibration object should be in the same plane as the organ to be measured.

Calibration precision directly depends on the length of the calibration distance: At longer calibration distances the measurement is more precise.

✓ (→ Page 97 *Calibrating images*)


- 1 Click the first point of calibration line in the image.
- 2 Click the second point of the calibration line in the image.

- 3 Enter the exact distance in mm.
- 4 Click **Calibration**.

### 4.8.6 Drawing arrows

To draw an arrow:

- 1 Click the ... icon and select **Arrow**.
- 2 Click the start point of the arrow in the image.
- 3 Click the end point of the arrow in the image.

(→ Page 97 *Deleting measurements*)

### 4.8.7 Storing an image

You can store a modified image and the original one is overridden.

If necessary, perform this work step: before saving the image, move the measurement, text item, or laterality label by dragging and dropping to a location where it does not cover up important details and where it is easy to read.

- ◆ Click this icon.

The image currently displayed is saved.

You can change the way to save the measurements you have performed.

(→ Page 122 *Image processing properties*)

### 4.8.8 Marking images

For sending images, you must mark them.

#### Marking the displayed image

- ◆ Click this icon.

Marked images are shown with a blue dot in the center at the left-hand side.

### Unmarking the displayed image

- ◆ Click this icon.

### Marking all images of the study

- ◆ Click this icon.

### Unmarking all images of the study

- ◆ Click this icon.

## 4.8.9 Displaying the exam protocol

All acquisitions performed on a patient are logged in the exam protocol.

- ◆ Click this icon.

| Exam Protocol   |                |     |    |           |        |                   |        |         |                     |             |      |                          |               | MULTIX Impact |           |
|-----------------|----------------|-----|----|-----------|--------|-------------------|--------|---------|---------------------|-------------|------|--------------------------|---------------|---------------|-----------|
| Patient Name    | 1811120812411  |     |    |           |        | Request ID        |        |         |                     |             |      | Number of Images         |               | 1             |           |
| Date of Birth   |                |     |    |           |        | Accession No.     |        |         | 20181120812411      |             |      |                          |               |               |           |
| Patient ID      | 20181120812411 |     |    |           |        | Study Description |        |         | Emergency           |             |      | Dose Area Prod.(uGy*cm2) |               | 0.00          |           |
| Patient Sex     | O              |     |    |           |        | Ref. Physician    |        |         |                     |             |      | Ref. Air Kerma(mGy)      |               | 0             |           |
| Ward            |                |     |    |           |        | Institution Name  |        |         |                     |             |      |                          |               |               |           |
| Irradiation no. | kV             | mAs | ms | Dose Lev. | MField | Cu.Filter         | Format | Density | Date and Time       | Application | DAP  | EX                       | Ref.Air Kerma | Focus         | Grid      |
| 1               | 40             | 0   | 0  |           | M      | NONE              |        | 0       | 2018-11-12 10:44:07 |             | 0.00 | 0.07                     | 0             | SMALL         | Universal |

## 4.8.10 Sending/archiving images

You can configure automatic sending/archiving of images.

DICOM properties

The images will be sent automatically to the configured sending/archiving nodes.

Only accepted images will be sent.

- ✓ At least one send/archive destination has been configured by Service.
- ✓ Images have been marked for sending. (→ Page 101 *Marking images*)
- 1 Select a **Target** from the list.
- 2 Click the **DICOM Dose Report** check box, if necessary.  
A DICOM Structured Report will be generated and sent with the image.
- 3 Click the **SC Image** check box, if necessary.

The exam protocol will be generated as a DICOM Secondary Capture image and this image will be sent with the image.

4 Click this icon.

The images are sent to the selected target in the background.

During transfer, an icon displayed bottom left, indicates the running transfer process.


Click this icon in the status bar to check the job status.

(→ Page 24 *Status of image transfers*)

## 4.9 Image hardcopy

In **Print mode**, you can view and arrange your images on virtual film sheets, adjust image sizes and view, select the hardcopy device/printer, and send the images to hardcopy.

- (→ Page 104 *Arranging film sheets*)
- (→ Page 106 *Changing the image size on the film sheet*)

- (→ Page 108 *Changing the image view on the film sheet*)
- (→ Page 108 *Printing images*)

You can configure automatic filming of images.

DICOM properties

The images will be sent automatically to the configured hardcopy camera.

Only accepted images will be filmed.

### 4.9.1 Arranging film sheets

In **Print** mode, your images are displayed on virtual film sheets. You can select the layout for the film sheets. You can adjust the size and view of each image.

#### Selecting and displaying a film sheet

- ◆ Click the film sheet.

#### Changing to portrait/landscape

- ◆ Click this icon.

#### Changing the film sheet layout

- 1 Click this icon.

- 2 Click the bottom right square of the required layout, e.g. for a 2x3 layout, click the square in the 2nd column and 3rd row.


### Creating a child layout

- 1 Click this icon.

- 2 Click the bottom right square of the required child layout, e.g. for a 2x2 child layout, click the square in the 2nd column and 2nd row.

Example

### Rejecting a child layout

- ◆ Click this icon.

### Repacking film sheets

- ◆ Click this icon.

Empty film sheets and segments between images are cleared.

All images are adjusted to the segment.

### Changing the order of images on the film sheets

- 1 Click this icon.
- 2 Drag and drop the images to the required positions.
- 3 Click the icon again.

### Adding a film sheet

- ◆ Click this icon.

### Deleting a film sheet

- ◆ Click this icon.

### Deleting an image from a film sheet

- 1 Click the image.
- 2 Click this icon

Click an image from the **Images** list, drag and drop it to the film sheet.  
The image is added to the film sheet.

#### 4.9.2 Changing the image size on the film sheet

##### CAUTION

Measurements made on an image that is not true-size

##### **Incorrect diagnosis possible**

- ◆ Do not assume that images are displayed in true-size on the monitor or film.
- ◆ Use the distance function for critical measurements.

**CAUTION**

Misinterpretation of the image size on the film

**Wrong measurements used for further treatment (e.g. prosthesis selection)**

- ◆ Check the scale shown on the film.
- ◆ Calculate any distances measured on the film to correspond with this scale.

Distortions are possible with projection images.

**Zooming an image continuously**

- 1 Click an image and then click this icon.
- 2 Drag the mouse cursor in the image up or down.

**Panning an image**

- 1 Click this icon.
- 2 Press the left mouse button and move the mouse cursor in the image.

**Adjusting to the segment**

- 1 Click an image.
- 2 Click this icon.

All the images on a film sheet are adjusted to the segment.

**Adjusting the same magnification**

- ◆ Click an image and then click this icon.

All the images on a film sheet are adjusted to the same magnification.

**Applying true size**

- ◆ Click this icon.

True size is applied to all the images on a film sheet.

### 4.9.3 Changing the image view on the film sheet

#### Rotating an image by $\pm 90$ degrees

- ◆ Click one of these icons.

#### Hiding/showing image text in preview layout

#### CAUTION

Patient name not displayed on film sheet

#### Incorrect diagnosis basis because of patient-image mismatch

- ◆ Be careful when configuring image text for the film sheets.

If the patient name is removed in the configuration, for example, to enable an anonymized film sheet, be sure to reset this setting immediately after use to make sure that diagnostic images are printed with the complete patient identification.

- ◆ Click this icon.  
The image texts of all the images on a film sheet are hidden or shown.

#### Inverting gray scales

- ◆ Click this icon.  
The gray scale of all the images on a film sheet are inverted.

### 4.9.4 Printing images

- ✓ **Print mode**
- ✓ At least one hardcopy device/printer has been configured by Service.
- ✓ Film sheets have been arranged, images have been adjusted in size and view.
- 1 Select a **Printer** from the list.
- 2 Select a **Film Size** from the list.
- 3 Enter the number of **Copies**.
- 4 Click the **Anonymous** check box, if necessary.

If checked, no patient information will be printed on the film sheets.

- 5 Click the **Print exam protocol** check box, if necessary.

The exam protocol will be generated as a DICOM Secondary Capture image and this image will be printed in a separate film sheet. (The layout for the exam protocol only can be 1x1.)

6 Click this icon.

The images are sent to the selected hardcopy device/printer in the background.

During transfer, an icon displayed bottom left, indicates the running transfer process.

Click this icon in the status bar to check the job status.

(→ Page 24 *Status of image transfers*)

## 4.10 Study/patient administration

- After examination or importing, the study/patient and image data is stored in the **Examined Patients** list.
  - (→ Page 26 *Study lists*)
  - (→ Page 26 *Accessing the study lists*)
  - (→ Page 28 *Searching for studies/patients*)
  - (→ Page 30 *Selecting a study/patient*)
  - (→ Page 109 *Displaying the images of studies*)
  - (→ Page 34 *Correcting study data*)
  - (→ Page 110 *Importing studies/images*)
  - (→ Page 112 *Exporting studies/images*)
  - (→ Page 114 *Sending/archiving studies/images*)
- After archiving the studies, you can delete them.
  - (→ Page 35 *Deleting studies*)
  - (→ Page 36 *Protecting/unprotecting studies*)
- If configured, studies are moved to the **Recycle Bin** when deleting them.
  - (→ Page 121 *General settings*)
  - (→ Page 116 *Managing studies in the Recycle Bin*)
- Information about DICOM archived studies is kept in the **Archiving Information** list.
  - (→ Page 116 *Retrieving studies/images*)

### 4.10.1 Displaying the images of studies

- ✓ **Patient mode**
  - ✓ **Pre-registered Patients** (with images acquired), **Examined Patients**, **Archiving Information**, or **Recycle Bin** (with images acquired) list
- 1 Click the study to display (to select multiple studies, use the Shift key or click the check boxes of the desired studies).

2 Click this icon.

### 4.10.2 Importing studies/images

You can import images in private format (\*.ecm) and DICOM (\*.dcm) CR or DX format (CR = Computer Radiography, DX = Digital X-ray).

Images in private format (\*.ecm) can only be imported together with DICOM format files (\*.dcm). If private format images exist on the storage medium, the private image will be imported with DICOM format images. The imported images can be processed in **Postprocessing**, **View**, and **Print** mode.

If private format images do not exist on the storage medium, only DICOM format images will be imported. The imported images can only be processed in **View** and **Print** mode.

An import is only possible when the same study does not exist in the current database.


- ✓ **Patient mode**
- ✓ **Examined Patients list**
- ◆ **When importing from CD/DVD:** Place a CD/DVD in the drive.
- or –
- When importing from USB memory:** Plug a USB memory device into the USB jack.

### Importing single images

- 1 Click this icon.
- 2 Browse to the drive/folder which contains your images.

- 3 Click the image file (\*.dcm).

(Images in private format (\*.ecm) are not displayed.)

– or –

Press the **Ctrl** or **Shift** key and click the image files (\*.dcm).

Start import by clicking **OK**.

Click **Abort** to cancel the import.

When the transfer is finished, a confirmation dialog is displayed. Click **OK** and the **Examined Patients** list is updated automatically.


Transfer can take some time depending on the volume of data.

### Importing a folder of images

- 1 Click this icon.

- 2 Browse to the drive/folder which contains your images.

- 3 Click the folder to import.

- 4 Click **OK**.
- 5 Start import by clicking **OK**.

A confirmation dialog is displayed if the folder contains more than 5 DICOM format files (\*.dcm).

Click **OK** to continue.

Click **Abort** to cancel the import.

When the transfer is finished, a confirmation dialog is displayed. Click **OK** and the **Examined Patients** list is updated automatically.

Transfer can take some time depending on the volume of data.

### 4.10.3 Exporting studies/images

You can export images in various formats:

- DICOM DCM
- DICOM DCM and ECOM ECM (private format)
- TIF
- BMP
- JPG

Patient information is only exported in DICOM format.

Measurement marks are only possible for export once they are burned into the image data.

(→ Page 122 *Image processing properties*)

During export, a DICOM viewer is automatically copied to the CD/DVD/USB.

Start the **syngo FastView** and select **Open File...** from the **Patient** menu to review the images stored in the CD/DVD/USB.

The exam protocol and DICOM Structured Report will be generated and exported automatically with the images.

- ✓ Patient mode
- ✓ Examined Patients list

### Exporting to the file system

- 1 When exporting to USB memory: Plug a USB memory device with enough free storage capacity into the USB jack.

Wait a moment.

- 2 Click the study to export in the **Examined Patients** list (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).


- 3 Click this icon.

To select the export folder:

- 4 Click **Browse**.
- 5 Browse to the drive/folder which shall contain your images.
- 6 Click the folder to export to.

- 7 Click **OK**.
- 8 Select the **File Format**.
- 9 Click **Create DICOM Dlr** for DCM and ECM format, if necessary.
- 10 Click **Anonymous**, if necessary.
- 11 Click **Export**.

You can click **Cancel** to abort the export.

- 12 When the export is finished, a confirmation dialog is displayed.

Click **OK**.

To view the files which are already stored:

- 13 Click **Open** and browse to the drive/folder which contains your images.

### Exporting to CD/DVD

- 1 Place a writable CD/DVD with enough free storage capacity in the drive.
- 2 Wait a moment.
- 3 Click the study to export in the **Examined Patients** list (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).
- 4 Click this icon.

- 5 Select the **File Format**.
- 6 Change other **Burn Settings**, if required.
- 7 Click **Burn**.

#### 4.10.4 Sending/archiving studies/images

You can configure automatic sending/archiving of images.

DICOM properties

The images will be sent automatically to the configured sending/archiving nodes.

Only accepted images will be sent.

**CAUTION**

Network connection may not be permanently available.

**Risk that images cannot be retrieved when they are needed.**

**Risk that images are lost.**

- ◆ Always retrieve patient and image data before starting examination.
- ◆ Always check transmission of images and restart transfer, if necessary.

**CAUTION**

The hard disk of the imaging system is not suitable for image archiving.

**Risk of data loss due to system failure (for example, defective hard disk)**

- ◆ Always archive patient and image data on suitable archive media as soon as possible.

- ✓ **Patient mode**
- ✓ **Examined Patients list**
- ✓ At least one send/archive destination has been configured by Service.
- 1 Click the study to send/archive (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).
- 2 Click this icon and select the image(s) you want to send/archive.

- 3 Click this icon.
- 4 Select a target or multiple targets from the list.
- 5 Click the **DICOM SR** check box, if necessary.  
A DICOM Structured Report will be generated and sent with the image.
- 6 Click the **SC Image** check box, if necessary.  
The exam protocol will be generated as a DICOM Secondary Capture image and this image will be sent as well.
- 7 Click **Send Image**.  
The images are sent to the selected target in the background.

### 4.10.5 Retrieving studies/images

Retrieve is only possible when the studies/images are no longer in the system.

#### CAUTION

Network connection may not be permanently available.

**Risk that images cannot be retrieved when they are needed.**

**Risk that images are lost.**

- ◆ Always retrieve patient and image data before starting examination.
- ◆ Always check transmission of images and restart transfer, if necessary.

- ✓ **Patient** mode
- ✓ **Archiving Information** list
- ✓ An archive has been configured by Service.
- 1 Select an archive from the **DICOM Node** list.

To test the connection:

- 2 Click this icon.
- 3 Define the search criteria, and click **Search**.

The studies meet the search criteria will display in the **Archiving Information** list.

(→ Page 28 *Searching for studies/patients*)

- 4 Click the study to retrieve in the **Archiving Information** list (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).

- 5 Click this icon.
- The images are imported in the background.

Retrieved images cannot be postprocessed anymore.

### 4.10.6 Managing studies in the Recycle Bin

When the **Recycle Bin** has been enabled, the deleted studies are moved to the **Recycle Bin**.

(→ Page 121 *General settings*)

The **Recycle Bin** list shows the deleted studies/images.

You can restore studies from the **Recycle Bin** or finally delete them.

Studies in the **Recycle Bin** are not deleted automatically.  
Make sure to empty the **Recycle Bin** regularly.

### Restoring studies from the Recycle Bin

- 1 Click the study to restore in the **Recycle Bin** list (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).

- 2 Click this icon.  
A confirmation dialog appears.
- 3 Click **Yes**.

### Deleting studies from the Recycle Bin

- 1 Click the study/studies to delete in the **Recycle Bin** list (to select multiple studies, use the **Shift** key or click the check boxes of the desired studies).

- 2 Click this icon.  
A confirmation dialog appears.
- 3 If you want to finally delete the selected studies: Click **Yes**.

### Emptying the Recycle Bin

- 1 Click this icon.  
A confirmation dialog appears.
- 2 If you want to finally delete all the studies in the **Recycle Bin**: Click **Yes**.

## 5 Administration of the imaging system

In **Service** mode, you can configure various settings, edit organ programs, perform statistical evaluations, and perform service functions.

- (→ Page 118 *User settings*)
- (→ Page 140 *Organ program and procedure editor*)
- (→ Page 158 *Service functions*)

### 5.1 User settings

In **Service** mode, you can configure various settings for your clinical practice.

- 1 Click one of the following buttons.
  - **Site Info.:** (→ Page 119 *Site information*)
  - **General Setting:** (→ Page 121 *General settings*)
  - **Processing Properties:** (→ Page 122 *Image processing properties*)
  - **Exam Properties:** (→ Page 124 *Examination properties*)
  - **Image Display:** (→ Page 125 *Image display properties*)
  - **Register Config.:** (→ Page 127 *Patient registration configuration*)
  - **DICOM Properties:**
  - **Disc Manager:**
  - **User Manager:** (→ Page 134 *User management*)
  - **Ortho setting:** (→ Page 132 *Settings for Ortho*)
- 2 Change the required settings.
- 3 Click **Apply** to accept your changes.
  - or –
  - Click **OK** to accept your changes and close the settings window.
  - or –
  - Click **Cancel** to reject your changes and close the settings window.

Some changes only take effect after system reboot.


---

### 5.1.1 Site information

- **Institution Name:** Name of your hospital or practice
- **Institutional Department Name:** Name of your department
- **Station Name:** Name of your ward
- **Address:** Address of your hospital or practice

### Login banner

#### Example


- **Display login banner:** Enable display of the login banner after startup
- **Title:** Title of the login banner
- **Content:** Content text of the login banner

### 5.1.2 General settings


- **Screen saver - Wait Time:** Set the time after which the screen saver will be activated. Possible settings: 0 ... 99 minutes, 0 = no screen saver
- **Use recycle bin:** Enable the Recycle Bin, i.e. studies are moved to the Recycle Bin when deleted.  
(→ Page 116 *Managing studies in the Recycle Bin*)
- **Use USB:** Enable the usage of USB memory devices.
- **Modality image:** Set the modality DICOM tag used for images.
  - **DX:** Digital X-ray modality
  - **CR:** Computer Radiography modality
- **Patient list:** Set single or multi study mode.
  - The **Single study mode** treats each study individually.
  - The **Multi-Study mode** treats all the studies of a patient together. When selecting a study, all studies of a patient are selected.  
(→ Page 26 *Study lists*)
- **The Default Position Of Mouse Cursor in Patient Mode:** Set the default position of the mouse cursor after entering Patient mode.
- **Sample of date and time format:** Set the format for date and time.
- **Show reject image:** Enable display of rejected images.
- **Save DICOM file with character set:** Set the language and encoding setting for DICOM files sent from the imaging system.

### Data in the study lists

Configuration of the displayed items, position, and sizes in the study lists

#### 5.1.3 Image processing properties


- **Show mask:** Enable the display of the detected ROI mask for image cropping.  
Image area outside the mask is displayed in black. The cropping mask is displayed with a white frame.
- **Show shutter:** Enable the display of the image shutter immediately after acquisition.  
When it is not checked, after acquisition no shutter is displayed. When manually cropping the image, the shutter will be displayed with a blue frame.

- **Fit to window when shuttered:** Enable automatic resize of image to the image display area when shuttered.
- **Measurement mark mode:** Select storage mode for measurement graphics in the image.
  - **Burn in image data:** measurement graphics are burned into the image data and will be archived with the image
  - **Image overlay:** measurement graphics will not be exported or printed with the images
- **Display image receptor size in:** Select which unit (cm or inch) to use for the image receptor format.
- **Tag font setting:** Select font and default font size for texts and labels in the image.
- **Distance shutter position:** Adjust the shutter position outward from the auto-cropping result, in pixels (0 to 200).  
0 means initializing shutter position by auto-cropping result.

### Configuring predefined texts/comments

You can save frequently used texts/comments as templates and select them from the list as needed.

- 1 To add a text item:  
Enter the text in the **Add Comment:** field and click **Add**.  
The text is added to the **Comments** list as a new entry.  
To delete a text item:
- 2 Select an item in the **Comments** list and click **Delete**.

## Configuring reject reasons

You can save frequently used reject reasons as templates and select them from the list as needed.


To add a reject reason:

- 1 Enter the text in the **Add Reason:** field and click **Add**.  
The text is added to the **Reject Reasons** list as a new entry.

To delete a reject reason:

- 2 Select an item in the **Reject Reasons** list and click **Delete**.

### 5.1.4 Examination properties

- **Show Predefined Demo Image:** Enable simulation of image acquisition without connection to the real X-ray system (for demonstrations, not for clinical use)

If enabled, an acquisition can be simulated by clicking the X-ray icon in **Examination** mode.

- **Heat Unit Test:** Enable the simulation of a high-dose X-ray shot

This test is required for the inspection of RAD systems in hospitals.

In **Examination** mode, press **Ctrl + H**.

The imaging system simulates a high dose X-ray shot and the heat unit display should be updated to **62%**. After restart of the imaging system, the test is disabled again.

- **Auto Accept Image:** Enable automatic acceptance of each acquired image.

(→ Page 55 *Releasing an exposure*)

- **Procedure Order:** Define the sort order for the procedure/organ program schedule:

- **Body Part:** Sort order for body parts
- **Protocol:** Sort order for procedures
- **View Position:** Sort order for organ programs of a procedure

Select an item and then click the up or down button to change the order.

### 5.1.5 Image display properties


#### CAUTION

Patient name not displayed on film sheet

**Incorrect diagnosis basis because of patient-image mismatch**

- ◆ Be careful when configuring image text for the film sheets.

If the patient name is removed in the configuration, for example, to enable an anonymized film sheet, be sure to reset this setting immediately after use to make sure that diagnostic images are printed with the complete patient identification.

- 1 Click a text item.
- 2 Enter/change the **Prefix**, **Postfix**, and **Default value**, if necessary.
- 3 Change the **Font Property Information**, if required.
- 4 Rearrange the text items via drag and drop, if necessary.


The patient name on the image is displayed in the name format from DICOM, i.e. Last Name^First Name^Middle Name.

- 5 To discard changes and reset to the defaults: Click **Factory Reset**.

**Label "left/right":**


- 6 Select text for laterality labels, depending on regional practice.
  - **L/R:**
    - Display **L** (English: Left or German: Links) for left label
    - Display **R** (English: Right or German: Rechts) for right label
  - **G/D:**
    - Display **G** (French: Gauche) for left label
    - Display **D** (French: Droite) for right label
  - **SIN/DX:**
    - Display **SIN** (Latin: sinister) for left label
    - Display **DX** (Latin: dexter) for right label

### 5.1.6 Patient registration configuration

- **Auto Generate PID:** Enable automatic generation of patient identification number during patient registration
- **Patient ID Auto Increment:** Enable automatic increment of patient identification number during patient registration
- **Current PID:** Start for automatic increment of patient identification number
- **Name-display format:** Define display of patient name components
- **Patient category:** Set the age for pediatric patients

#### Input fields configuration


**Input Fields:** Define the input field available during patient registration

- **Visible:** Item will be shown.
- **Mandatory:** Item must be filled out. (Item will be indicated by an asterisks \*.)

(→ Page 31 *Manual registration*)

### 5.1.7 DICOM properties

- **Worklist query:**
  - **Waiting time:** Update interval for the **Pre-registered Patients** list.
  - **Specific Character Set:** Selection of specific character set.
  - **Use Unicode (UTF-8/ISO\_IR 192) for DICOM character encoding:** Enable unicode for DICOM character encoding.
- **Disable auto RIS query when entering patient mode**
  - **Yes:** The **Pre-registered Patients** list will not be updated automatically when the system switches to the **Patient** mode.
- **Close study method:**
  - **Auto on entering patient mode:** Enable automatic MPPS close when switching back to the **Patient** mode.
  - **Manually with dialog:** Enable manual MPPS close when switching back to the **Patient** mode.

- **Storage commitment:**
  - **Prepare Commitment:** The time the imaging system is waiting after a study / images have been sent to an archive server with storage commitment.
  - **Commitment timeOut:** Time-out for storage commitment
  - **Number of attempts:** The imaging system performs the selected number of retries when the archive server is not available.


## Autosending configuration

- **Auto-sending RAD Images:** Selection if and when automatic sending is performed:
  - **Off:** No images will be sent automatically
  - **When study is closed:** Images will be sent to the autosending destination after study close.
- **Default autosending destination:** Configure up to three hosts which shall be used as the default send destinations.
- **Autosending dose information:**
  - **SC Image:** The exam protocol will be autosent as a secondary capture image.
  - **DICOM Dose SR:** The exam protocol will be autosent as a DICOM structured report.
  - **Destination:** Select destination for autosending dose information.

It is also used for manual sending, when predefined send destinations are used.

## Autofilming configuration

- **Autofilming RAD Images:** Selection if and when automatic filming is performed:
  - **Off:** No images will be sent to the autofilming hardcopy camera automatically
  - **When study is closed:** Images will be sent to the autofilming hardcopy camera after study close.
- **Default Autofilm Camera:** Configure the default hardcopy camera which shall be used as the default hardcopy camera.
- **Background color for exam protocol:** Set the background color for the exam protocol DICOM Secondary Capture image.
- **Magnification type:** Set the interpolation type by which the printer magnifies or decimates the image in order to fit the image in the image box on a film.

### 5.1.8 Disk management

#### Auto Delete:

- **Enable auto delete:** Activate the automatic image deletion function  
If enabled, you can define the rules for automatic deletion.
- **Free space after deleting more than:** Minimum free space on the hard disc for activation of automatic deletion
- **Activation time:** Activation time for automatic deletion (hh:mm:ss)
- **Activation date:** Activation day for automatic deletion, **Daily** or weekly with selection item for day
- **Conditions:** Condition for automatic deletion of images (arbitrary combination of following conditions)
  - **Image sent:** Delete all images which have been sent to a configured archive node (with or without storage commitment).
  - **Image filmed:** Delete all images which have been filmed.
  - **Image stored on CD/DVD:** Delete all images which have been stored on CD/DVD or USB memory.

In order to enable automatic deletion, at least one condition must be set.

**Disk space monitoring:** Display a warning message when a disk fill level is reached

- **Disk monitoring interval:** Check the disk fill level every x minutes
- **Warning if free space is less than:** Disk fill level for warning message

Values can be set  $\geq 3$  G.

When the free disk space reaches the warning level, an automatic image deletion confirmation dialog is displayed.

- Click **OK** to start the deletion.
- Click **Cancel** to continue the examinations.

- **Locked if free space is less than:** Disk fill level for blocking acquisition

Values can be set  $\geq 2$  G.

When the free disk space reaches the lowest limit level, acquisitions are blocked and an automatic image deletion is required.

Click **OK** to start the deletion.

You can also clear the disk manually.

(→ Page 158 *Clearing study and image data*)

## 5.1.9 Settings for Ortho

### General setting


**Distance shutter position:** Adjust the shutter position inward or outward from the auto-cropping result, in pixels.

0 means initializing shutter position by auto-cropping result.

## Organ detection settings

You can adjust the collimation width factor and offset for Ortho acquisitions with auto full spine and long leg collimation.

**Full spine AP/PA width factor:** Collimation ratio (30% to 150%)

**Full spine upper offset:** Offset on the upper side for auto full spine (-200 mm to 200 mm)

**Full spine lower offset:** Offset on the lower side for auto full spine (-200 mm to 200 mm)

**Long leg upper offset:** Offset on the upper side for auto long leg (-200 mm to 200 mm)

**Long leg lower offset:** Offset on the lower side for auto long leg (-200 mm to 200 mm)

(→ Page 134 *Setting instruction*)

## Setting instruction


Illustration for offset

## 5.2 User management

The imaging system handles different users and is able to identify and distinguish between every specific user.

(→ Page 134 *Changing the password*)

### 5.2.1 User roles and privileges

The following roles are available in the imaging system as default.

| User role         | Privileges                                                     |
|-------------------|----------------------------------------------------------------|
| Tech              | Basic examination, post-processing and documentation operation |
| Emergency         | Limited examination operation                                  |
| Mgr               | Clinical administration                                        |
| Siemens           | Siemens Healthineers service                                   |
| IS_Administrators | Local user administration                                      |
| IS_Auditors       | Access to audit trail files                                    |

You can change the roles as required.


(→ Page 135 *Setting up users, roles, and permissions*)

### 5.2.2 Changing the password

✓ Service mode

- 1 Click **User Manager**.

- 2 Enter your old password.
- 3 Enter your new password.
- 4 Re-enter your new password.
- 5 Click **Save**.

### 5.2.3 Setting up users, roles, and permissions

- ✓ You are logged in with a user account which has the permission to access these functions.
- ✓ **Service mode**

#### CAUTION

Changing users, groups, and settings


**By using the administration tools, it is possible to modify the imaging system so that acquisition is no longer possible.**

- ◆ It is the responsibility of the administrator to set up users and groups correctly.
- ◆ Observe the notes in this chapter.

- ◆ Click **User Manager**.

### Changing user accounts

- 1 Click **User**.

To add a new user account:

- 2 Click **New** and enter a **User Name**.

Enter a **Password**.

Assign a **Role Name**.

To change a user account:

- 3 Click a user account in the **User List** and change the **User Name**.


Assign a **Role Name**, if necessary.

To delete a user account:

- 4 Click a user account in the **User List** and click **Delete**.
- 5 Click **Save**.

### Assigning permissions to roles

- 1 Click **Role**.

To add a new role:

- 2 Click **New** and enter a **Role Name**.

Click the permissions to assign in the **Permission** list.

To change a role:

- 3 Click a role in the **Role Name** list.

Click the permissions to assign in the **Permission** list, if necessary.

To delete a role:

- 4 Click a user account in the **Role Name** list and click **Delete**.

A confirmation dialog is displayed.

- 5 Click **Yes**.


Click **Save** to confirm the changes.

## Setting up security options

You can reset the password for any user.

- 1 Click **Cyber Security**.

- 2 Select a **User Name**, enter a new **Password**, and click **Reset**.

For the first login after password reset, the user will be asked to change the password.

To set up password rules:

- 3 Enter the criteria in the **Password Rules** list.

### Setting up Remote Interface (option)

You can register or authenticate the Remote Interface.

For more information, refer to the Remote Interface Operator Manual.

#### Registering the Remote Interface

(You are logged in with the administrator role.)


- 1 Click **Remote Interface**.

- 2 Click **New**.
  - 3 Enter the **Mobile name** and **Device UID**.
  - 4 Click **Save**.
- The mobile device is listed in the **Mobile list**.

#### **Deleting the Remote Interface** (You are logged in with the administrator role.)

- 1 Click **Remote Interface**.

- 2 Select the **Remote Interface** from the **Mobile list** and click **Delete**.  
A message appears asking if you want to delete the **Remote Interface**.
  - 3 Click **Yes** to confirm the deletion.
- The **Remote Interface** is removed from the **Mobile list**.

### Authenticating the Remote Interface

(You are logged in with one of the user roles: Tech, Mgr or Siemens.)

#### 1 Click Remote Interface.

- 2 Scan the QR code with your Remote Interface  
Click Refresh to update the QR code if necessary.

## 5.3 Organ program and procedure editor

The organ program and procedure editor, referred to as the **OGP Editor**, allows to define organ programs (OGPs), define procedures that consist of one or more organ programs, and select parameters for filming/printing. Additionally, administrative tasks can be performed in the OGP Editor.

(→ Page 143 *The OGP Editor*)


### 5.3.1 General information

#### Access modes/access control

Following modes for accessing the OGP Editor are available:

- In **Customer mode**, you have access rights for the customer data, and additionally for viewing the default procedures and organ programs. Access to quality control and detector calibration organ programs is not possible.
- In **Mgr mode**, you have access for viewing and modifying the default procedures and organ programs.
- In **Siemens mode**, special Siemens Healthineers image quality experts have access for viewing and modifying the default procedures and organ programs as well as quality control and detector calibration organ programs.

To open the OGP Editor you are asked to enter a password. The customer mode password is site specific, you can set and modify it.

(→ Page 142 Opening the OGP Editor)

## Online/offline mode

The OGP Editor is a tool for viewing/editing of organ programs and procedures of the OGP database, it runs on the imaging system.

Offline editing of the OGP database is also supported. The OGP Editor can run on a separate Windows PC without the main imaging system application.

## OGP database

The OGP Editor handles a database, referred as the OGP database, containing the definitions of the procedures and programs.

## Import/export

It is possible to import a OGP database into the OGP Editor via any local or external device, e.g. local hard disk, CD/DVD, USB memory. Each import replaces the current database and cannot be undone.

An import is only possible for a database of the same system type. Depending on licence settings, some procedures/OGPs may become inactive after transfer. Inactive procedures/OGPs cannot be selected.

It is also possible to export the database to USB, CD, DVD, and local hard disk for external modification or transfer to other systems or to save a copy of the current database (backup).

## Automatic backup

Backups of the OGP database are automatically performed. The automatic backup runs in the background after each editing session, when closing the OGP Editor.

Up to 20 backups are kept in a folder on the hard disk.

A restore of a database stored via automatic backup can be done with the import function.

## Organ programs and procedures

The names of the procedures and organ programs depend on the set language of the user interface:

- Siemens Healthineers default procedures and organ programs are language dependent.

In other words, each organ program or procedure has a different name for each configurable language, but the parameters are the same for all languages.


- Each procedure can contain up to 50 organ programs.

### 5.3.2 Opening the OGP Editor

- ✓ **Service mode**

- 1 Click **OGP Editor**.

The password dialog is displayed.


- 2 Enter or select an OGP Editor account from the list.
- 3 Enter the correct password.

The factory-set password for customer mode is: 137408

- 4 Click **Login** or press the **Enter** key.

### 5.3.3 The OGP Editor

#### Example

The **OGP Editor** presents the following tabs:

- **Procedure:** Define and edit procedures.  
(→ Page 144 *Procedure data*)  
(→ Page 146 *Assigning organ programs to procedures*)
- **Organ Program:** Define and edit organ programs and parameter subsets.  
(→ Page 149 *Organ program parameters*)
- **Configuration:** Define and edit system configuration data, e.g. language configuration, password setting. Functions for backup/restore of the OGP database are also located here.  
(→ Page 155 *Configuration parameters*)  
(→ Page 143 *Editing procedures*)

Depending on your system and the installed options not all parameters may be displayed or applicable.

### 5.3.4 Editing procedures

- ✓ OGP Editor in Procedure tab

#### Creating a new procedure

- 1 Click **New**.

- 2 Enter/select the required data in the fields.
- 3 Select at least one related organ program.
- 4 Click **Save**.

### Deleting a procedure

- 1 Select the procedure to delete.
- 2 Click **Delete**.  
A confirmation dialog appears.
- 3 Click **Yes**.


### Adopting changed settings

- ◆ Click **Save**.

### Rejecting changed settings

- ◆ Click **Cancel**.

## 5.3.5 Procedure data

Example

In the **Procedure** tab, the **OGP Editor** displays all available procedures and organ programs:

- In the **Procedure list**, all procedures are listed.  
The order in a body part is based on the settings in the **Examination Properties** in the user settings.  
Newly added procedures are listed at the bottom.
- The **Procedure OGP list** shows the organ programs which are assigned to the selected procedure.
- In the **OGP list**, all organ programs are listed.  
You can double click an organ program in the **OGP list** to edit it.

The procedures and organ programs in the lists can be filtered by the **Group**, i.e. body part.

The **Procedure information** displays assigned data, you can enter/change them:

- **Code:** Code of the procedure  
It is the identifier of a procedure. Each procedure has its own code. Do not use same code for different procedures!
- **Group:** Body part for the procedure
- **Print style:** Default layout for filming/printing  
Click **Details** to display the **Print style details**.

- **Film size:** Film sheet size (length x width)
- **True size:** Zoom image to true size yes/no
- **Film orientation:** Orientation of film sheet portrait/landscape
- **Film format:** Format/layout of images on film sheet
- **Often used:** Mark as a favorite procedure.  
The procedure marked as **Yes** will be displayed in **Often used** category in **Procedure Information**.  
(→ Page 32 *Entering procedure information*)
- **Name:** Name of the procedure
- **Category:** Patient category for the procedure - **Adult** or **Pediatric**

- **Description:** Description for the procedure
- **Enable:** Mark as an active procedure.

### Assigning organ programs to procedures

- ◆ Select a procedure in the **Procedure list**.  
The assigned organ programs are shown in the **Procedure OGP list**.

- To add an organ program**
- 1 Select an organ program in the **OGP list** and click <.
  - 2 Click **Save**.

For Ortho procedures, only the organ programs that belong to the same group can be used.

- To remove an organ program**
- 1 Select an organ program in the **Procedure OGP list** and click >.
  - 2 Click **Save**.

### 5.3.6 Editing organ programs

- ✓ **OGP Editor in Organ Program tab**

#### Creating a new organ program

- 1 Click **New**.
- 2 Enter/select the required data in the fields.
- 3 Enter/select the required parameters.
- 4 Link one to four X-ray parameter templates (according to your system configuration) and one image processing template to the organ program.  
Select an **Acq type** as default workplace.
- 5 Click **Save**.

An Ortho procedure consists of up to 4 Ortho organ programs.

After saving an Ortho organ program successfully, the following message appears:  
<"Do you want to create another Ortho OGP?">

Click **Yes** to create another Ortho organ program in the same group.

#### Deleting an organ program

- 1 Select the organ program to delete.  
Make sure that the organ program is not assigned to any procedures.  
(→ Page 146 *To remove an organ program*)

- 2 Click Delete.

A confirmation dialog appears.

- 3 Click Yes.

### Adopting changed parameters

- ◆ Click Save.

### Rejecting changed parameters

- ◆ Click Cancel.

### Printing the organ program parameters

- ✓ A paper printer has been configured.
- ◆ Click Print.

### Managing X-ray parameter templates

- ✓ OGP Editor in Organ Program tab

#### Removing a X-ray parameter template


- 1 Select the Acq type (workplace template) to remove from the organ program.
- 2 Click Remove.  
Select an Acq type as default workplace.
- 3 Click Save.

#### Assigning a new X-ray parameter template to an organ program

- 1 Click Edit to open the X-Ray Parameter template tab.
- 2 Select the required Acq type (workplace template) from the template list.
- 3 Click Apply.

The OGP Editor dialog is displayed.

A confirmation dialog appears.

- 4 Click OK.
- 5 Enter/select the required data in the fields.

- 6 Click **Save**.

A confirmation dialog appears.

- 7 Click **OK**.

#### Creating a new X-ray parameter template

- 1 Click **Edit** to open the **X-Ray Parameter** template tab.
- 2 Click **New**.
- 3 Enter/select the required data in the fields.
- 4 Click **Save**.

A confirmation dialog appears.

- 5 Click **OK**.

Newly added templates are listed at the bottom of the **template list**.

#### Deleting an X-ray parameter template

- 1 Click **Edit** to open the **X-Ray Parameter** template tab.
- 2 Select the required **Acq type** (workplace template) from the **template list**.
- 3 Click **Delete**.

A confirmation dialog appears.

- 4 Click **Yes**.

A confirmation dialog appears.

- 5 Click **OK**.

If the template is used by some organ programs, deletion is not possible.  
Remove the template from the organ programs listed in the **Related OGP** list.

## Managing image processing parameter templates

### ✓ OGP Editor in Organ Program mode

#### Assigning a new image processing parameter template to an organ program

- 1 Click **Edit** to open the **Image processing parameter** template tab.
- 2 Select the required template from the **template list**.
- 3 Click **Apply**.

The **OGP Editor** dialog is displayed.

- 4 Select the default flavor and click **Save**.

#### Adding a new image processing parameter template

- 1 Click **Edit** to open the **Image processing parameter** template tab.
- 2 Click **New**.
- 3 Enter/select the required data in the fields.
- 4 Click **Save**.

A confirmation dialog appears.

- 5 Click **OK**.

Newly added templates are listed at the bottom of the **template list**.


#### Deleting an image processing parameter template

- 1 Click **Edit** to open the **Image processing parameter template** tab.
- 2 Select the required template from the **template list**.
- 3 Click **Delete**.  
A confirmation dialog appears.
- 4 Click **Yes**.  
A confirmation dialog appears.
- 5 Click **OK**.

If the template is used by some organ programs, deletion is not possible.  
Remove the template from the organ programs listed in the **Related OGP** list.

### 5.3.7 Organ program parameters

#### Example

- **Search:** Selection of organ program
  - **Body part:** Selection of body part to be examined
  - **OGP list:** Selection of organ program
  - **Show previous or next organ program**
- **Attribute:** Identification of organ program and miscellaneous parameters  
(→ Page 153 *Miscellaneous parameters*)

- **X-Ray Parameter:** Identification of X-ray parameter templates and their parameters  
(→ Page 150 *X-ray parameters*)
- **Image processing parameter:** Identification of image processing parameter templates and their parameters  
(→ Page 152 *Image processing parameters*)

## X-ray parameters

- **Acq type:** Workplace to be used.  
Each organ program may have more than one workplace templates.  
(→ Page 147 *Managing X-ray parameter templates*)
- **Default:** Default workplace for the organ programs  
You can select the default workplace for the organ program by clicking the corresponding check box and **Save** button.
- **Focus:** Select the focus setting of the X-ray tube: small/large

Other parameters can be changed in the **X-Ray Parameter** tab.

Click **Edit** to enter the template in order to change the parameters if necessary.

### Example

**Template list:** List of X-ray parameter templates

- **Name:** Name filter for parameter templates  
Enter a name in order to search in the templates.  
To select a template: Click an entry in the list.

- **Acq type:** Acquisition type filter for parameter templates
- **Related OGP:** List of organ programs using the parameter template

**Common attribute:**

- **ID:** Template identifier
- **Name:** Template name
- **Acq type:** Workplace type
- **SID:** Source to image receptor distance
- **Grid type:** Use of a grid and specification of grid type: **0, 115, Universal, 180**  
0 means that no grid is used.
- **Filter type:** Thickness of the copper filter: **0.0 mm, 0.1 mm, 0.2 mm, 0.3 mm**
- **Collimator center:** Alignment of collimation: **Top alignment, Center**
- **Collimator width:** Collimation width
- **Collimator length:** Collimation height
- **Direction:** Direction of Ortho acquisition: **Feet to Head or Head to Feet**
- **Tube height:** Tube height in cm defined for SmartPositioning function
- **Tube angle:** Tube angle defined for SmartPositioning function
- **Alignment:** Collimation alignment: **Center alignment, Top alignment**
- **SID tracking:** SID tracking function in table mode or wall mode (vertical projection):  
**Not active, Active, Lost**  
(only effective when the SmartPositioning function is configured)
- **Tracking:** Tube tracking function in table or wall mode: **Not active, Active, Lost**  
(only effective when the SmartPositioning function is configured)

**Condition attribute** - depends on selection of exposure mode:

- **Small/Medium/Large:** Patient size
- **Exp mode:** Selection of exposure mode
  - **1 pt:** 1-point technique (with Automatic Exposure Control, AEC):  
(not with free exposure)  
An acquisition is taken with operator specified kV value.  
**1 pt with ext. exp.:** 1-point technique with extended exposure (AEC with current reduction)
  - **2 pt:** 2-point technique:  
An acquisition is taken with operator specified kV and mAs values.
  - **3 pt:** 3-point technique:  
An acquisition is taken with operator specified kV, mAs and ms values.
- **kV:** Tube voltage [kilo volts]

- **ms:** Pulse width [milli seconds]
- **mAs:** Product of tube current and time [milli ampere seconds]
- **AEC field:** Selection of measuring fields (dominants) for organ adapted control of exposure on units.
- **Dose:** Exposure dose in  $\mu\text{Gy}$ .
- **Density:** To obtain darker or brighter images, the dose can be modified in steps of exposure points.
- **Tube load:** Maximum load for X-ray tube: 80%/100%
- **TOD:** Table-object distance


### Image processing parameters

- **ID:** Display of identifier of selected image processing template
- **Name:** Display of name of selected image processing template
- **Flavors:**
  - **Name:** Name of flavor in this template
  - **Default:** Default flavor

You can select the default flavor for the organ program by clicking the corresponding check box and **Save** button.

- **Parameters:** Display of image processing parameters

Click **Edit** to enter template subtab and change settings if necessary..

Example

- **Template list:** List of image processing templates
  - **Name:** Name filter for image processing templates  
Enter a name in order to search in the image processing templates.  
To select an image processing template: Click an entry in the list.
- **Related OGP:** List of organ programs using the processing template
- **Attribute:**
  - **ID:** Identifier of image processing template
  - **Name:** Name of image processing template
  - **Acq mode:** Acquisition mode
    - DR RAD
    - DR Cassette
    - DR RAD Ortho
- **Flavor setting:**
  - **Child flavors:** List of flavors in this image processing parameter template  
Child flavors display in the **Flavor Setting** tab in **Postprocessing** mode.
  - **Flavors:** List of available image processing flavors  
To add or remove a flavor, click the arrows.  
Click **Save** to adopt the change.
- **Parameters:** Display of image processing parameters


## Miscellaneous parameters

### Attribute:

- **ID:** Identifier of organ program
- **Name:** Name of organ program
- **Description:** Description for organ program
- **Acq mode:** Selection of acquisition mode
- **Body part:** Assignment to body part
- **Anatomic region:** Selection of anatomic region
- **Image laterality:** Laterality information of image content
  - **Not defined:** Image laterality not defined
  - **R** Right one of a paired organ can be seen in the image
  - **L:** Left one of a paired organ can be seen in the image
  - **U:** Unpaired organ can be seen in the image
  - **B:** Both paired organs can be seen in the image

- **Image flip and rotation:** Presetting for image flip or flip with rotation and the default position of flip indicator
- **Image rotation:** Presetting for image rotation in steps of 90 degrees
- **Image label:** Presetting for left/right label and default position of label
- **View position:** Position of the patient during examination:
  - **Not defined:** View position is not set
  - **AP:** Anterior - Posterior view
  - **PA:** Posterior - Anterior view
  - **RL:** Right Lateral view
  - **LL:** Left Lateral view
  - **RLD:** Right Lateral Decubitus view
  - **LLD:** Left Lateral Decubitus view
  - **RLO:** Right Lateral Oblique view
  - **LLO:** Left Lateral Oblique view
- **Target EXI:** Typical clinical exposure index value  
(→ Page 56 *Exposure index*)
- **View icon:** Selection of positioning graphics
- **Position guide:** Selection of positioning guide image
- **Demo Image:** Selection of demonstration image
- **Collimator no change:** Change collimation size:  
Click **Yes** to keep the collimator size of the last acquisition. Click **No** to use the collimator size set for this OGP
- **Orientation column/Orientation row:** Display of orientation marks in the image
  - **Not defined:** Orientation not defined
  - **A:** Anterior
  - **P:** Posterior
  - **L:** Left
  - **R:** Right
  - **H:** Head
  - **F:** Feet

### 5.3.8 Configuration parameters

#### Example

- **Language:** Language of the **OGP Editor** user interface.
- **Unit of length:** Unit for length measurements, for collimation size, and film size
- **Password:** (→ Page 155 *Changing the OGP Editor password*)
- **Software Information:** Display of OGP Editor and database versions and last updated date
- **Database management:** (→ Page 156 *Database management*)
- **Film size config:** Default film width/height  
To delete a film size, select the film width and height and click **Delete**.  
To add a film size, enter the film width and height and click **Save**.
- **Default acquisition type config:** Default acquisition type
  - **Original Acq Type:** Original default workplace
  - **Target Acq Type:** Desired default workplace
 Click **Switch** to change the default workplace for all organ programs.


The change of default workplace applies to all organ programs.

### Changing the OGP Editor password

- ✓ OGP Editor, Configuration mode

- 1 Click **Change password**.

- 2 Enter your old password.
- 3 Enter your new password.
- 4 Re-enter your new password.
- 5 Click **Save**.

## Database management

It is possible to import or export the OGP database via any local or external device (e.g. local hard disk or USB memory).

The use of OGP databases from products other than MULTIX Impact C is not supported and can result in insufficient image quality or loss of image.

OGP databases from other systems with a different software version are also not always compatible.

In case of doubt, contact Siemens Healthineers Service.

### ✓ OGP Editor, Configuration mode

#### Backing up the OGP database

- 1 Click **Backup**.  
A confirmation dialog is displayed.
- 2 Click **OK**.

Perform backup before you make major changes to the OGP database.

#### Exporting the OGP database

- 1 Click **Export**.
- 2 Select the destination drive and folder.
- 3 Click **Save**.  
A confirmation dialog is displayed.
- 4 Click **OK**.

### Importing an OGP database

Each import replaces the current database and cannot be undone.

- 1 Plug-in the USB memory, if necessary.  
Wait a moment.  
To open the file folder:
- 2 Click **Import** to replace the current database.
- 3 Select the source drive and folder.  
A list of all restorable databases is shown.
- 4 Click **Open**.  
The following message appears: <"Do you want to reserve local RIS mapping result?">
- 5 Click **Yes** to preserve your local RIS mapping result.  
– or –  
Click **No** to replace your local RIS mapping result by the imported database.  
A confirmation dialog is displayed.
- 6 Click **OK**.

### Merging an OGP database

Only organ programs not existing in the current database can be merged.

- 1 Plug-in the USB memory, if necessary.  
Wait a moment.  
To open the file folder:
- 2 Click **Merge** to merge the stored database with the current database.
- 3 Select the source drive and folder.  
A list of all restorable databases is shown.
- 4 Click **Open**.  
A confirmation dialog is displayed.
- 5 Click **OK**.

### 5.3.9 Closing the OGP Editor

- 1 Click the **X** button in the OGP Editor dialog.  
A confirmation dialog is displayed.
- 2 Click **Yes** to close.

Another confirmation dialog is displayed if changes are not saved. Click **Yes** to close without saving.

- or -

Click **No** to return to the dialog.

## 5.4 Service functions

Tool area:

- **OGP Editor:** (→ Page 140 *Organ program and procedure editor*)
- **Clean Archived Image:** (→ Page 158 *Clearing study and image data*)
- **FD Calibration:** (→ Page 159 *Detector calibration*)
- **FD Reconnect:** (→ Page 163 *FD reconnect*)
- **RIS Mapping:** (→ Page 163 *RIS mapping*)
- **Statistics:** (→ Page 168 *Statistics*)
- **Audit Trail:** (→ Page 169 *Audit trail*)
- **Test Image:** (→ Page 170 *Test images*)
- **Tube Load:** (→ Page 172 *Tube load*)

Service area:

- **PMS (Service Software):** (→ Page 173 *Service login*)
- **Remote Assistance:** (→ Page 173 *Remote assistance*)
- **Remote Access:** (→ Page 174 *Remote service access*)
- **DICOM Node:** (→ Page 175 *DICOM nodes*)
- **FD Manager:** for Siemens Healthineers service use only
- **Sensitivity:** (→ Page 180 *Sensitivity calibration*)

See also:

- (→ Page 23 *System messages*)

Additional information:

- (→ Page 172 *Providing an image for service*)
- (→ Page 172 *Providing log files for service*)
- (→ Page 181 *QR codes*)
- (→ Page 181 *HIPAA*)

### 5.4.1 Clearing study and image data

You can delete archived studies/images of a data range.

Studies and images can be archived via network, CD/DVD/USB memory or hardcopy.

✓ **Service mode**

1 Click **Clean archived image**.

2 Enter the **Start Date** for the studies to delete.

3 Enter the **End Date** for the studies to delete.

The number of studies is indicated.

4 If you want to delete all the studies of the selected date range now: Click **Begin Clean**.

Studies will be deleted, if a study fulfills the following conditions:

- At least one image has been archived via network, CD/DVD/USB memory, or hardcopy
- AND Study is not protected.
- AND Study is closed.

#DROC Disk management#

### 5.4.2 Detector calibration

Detector calibration is for adjusting the X-ray sensitivity of the flat detector.

(→ Page 160 *Performing FD calibration*)

#### Calibration reasons

A detector calibration must be performed regularly and:

- When the system is started for the first time at the system startup by a service technician.
- After detector has been exchanged.

- If exposure conditions such as distance between the detector unit and X-ray tube have been changed significantly.
- If exposure conditions have changed significantly.
- If a calibration warning appears.

The detector calibration should be done by a technician.

Acquisition will not be disabled if calibration is required. To acquire without calibrating is on your own risk. Image quality might not be sufficient.

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## Calibration requirements

Temperature in the room when performing calibration should be the same as in an ordinary exposure. Otherwise, you cannot obtain a uniform image.

Calibration can only be performed if the detector warm-up time of 10 minutes has elapsed.

Calibration must be done by exposing the whole detector area. Although calibration can be performed even by irradiating the detector unit partially, you will not be able to obtain a normal image for the part where X-ray has not been irradiated.

Do not put any lead mark on the detector unit, or place any object between the detector unit and the X-ray tube. Otherwise, object will appear on all the pictures.

## Performing FD calibration

The detector should be thermalized before calibration.

Make sure that there are no additional objects in the beam path, e.g. a filter.

- ✓ You are logged in with a user account which has the **FD Calibration** permission.
- ✓ **Service mode**

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- 1 Move the tube stand out of the tabletop and turn it to the floor.
- 2 Put the detector on the floor and turn it to landscape position.

Put the detector on a plastic foam sheet. The floor may be too cold.

- 3 Set the tube to maximum SID.  
At least 150 cm, 180 cm is recommended.
- 4 Measure the SID.
- 5 Center the detector in beam field and collimate the beam field to fully cover the active detector area.
- 6 Insert the 0.6 mm Cu filter.
- 7 Click **FD Calibration**.

**FD Calibration**

Calibration conditions

Detector: Venu1717XS\_B Serial number: VS36000211105190002

**Check FD parameters for calibration**

Instructions: Please ensure that no additional objects are in the beam path (e.g. filter)! Please wait while the system verifies detector parameters.

| Parameter            | Ready | Current value | Remarks                                                                                                                                                                                        |
|----------------------|-------|---------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| SID                  | NO    |               | Center the detector in beam field and collimate the beam field to fully cover the active detector area. Please set the tube to maximum SID (at least 130 cm) and enter the measured value[cm]. |
| Grid                 | n.a.  |               | Please remove the grid!                                                                                                                                                                        |
| Temperature constant | NO    | 36.98°C       | Temperature range [10,54]°C, Please wait 10min                                                                                                                                                 |
| Battery Capacity     | YES   | n.a.          | Battery capacity bigger than 50%                                                                                                                                                               |

Cancel OK

Please wait while the system verifies detector parameters.  
Warm up time may be needed to wait.


- 8 Enter the measured SID in the **Current value** box and click **OK**.

The following message appears: <"Is the following filter inserted: 0.6 mm Cu? Click 'Cancel' to abort the calibration">

- 9 Ensure that the 0.6 mm Cu filter has been inserted and click OK

The following message appears: <"Please check calibration switch ON.">

- 10 Turn on the calibration switch and click OK to continue.

The calibration switch can be found at the back side of the generator ON/OFF console.

- 11 Choose the calibration mode.

To perform the detector calibration in radiography mode, do not click the check box of **Long time exposure**.

– or –

To perform the detector calibration in long time exposure mode, click the check box of **Long time exposure** (not available in China market).

- 12 Ensure that nobody is in the examination room and click **Begin** to start automatic X-ray release.

The calibration is started.

Do not turn off the calibration switch when the calibration is in progress!

- 13 Wait until the calibration has finished.

The following message appears: <"The calibration has finished successfully. Click 'OK' to save the calibration.">

- 14 Click OK to save the calibration results.

- 15 Wait until the calibration results have been saved successfully.

The following message appears: <"Please close service switch, then restart system.">

- 16 Click OK.

- 17 Turn off the calibration switch and reboot your system.

It is very important to turn off the calibration switch for normal operation, otherwise, no radiation is possible!

If the calibration has not succeeded, the previous calibration remains valid.

### 5.4.3 FD reconnect

You can reconnect a disconnected detector, for example, after loss of WiFi signal or after recharging.

- ✓ Service mode
- Click FD reconnect.

The system starts connection to the detector.

### 5.4.4 RIS mapping

The imaging system provides a learning mode in order to map scheduled procedures of the RIS (radiological information system) with local procedure names.

Multiple RIS names can be mapped to one local procedure name.

One RIS name can also be mapped to multiple local procedure names.

The RIS mapping tool can be activated and inactivated. Access to this tool is possible either under password protection or without password.

The system requires one of three DICOM tags from RIS to perform RIS mapping:

- Requested procedure description (0032,1060)
- Scheduled procedure step description (0040,0007)
- Scheduled protocol code value (0040,0008)>(0008,0100).

Select one of the three tags for RIS mapping on the system.

All previous successful mapping results will be lost once the DICOM tag has been changed.


After successful mapping, when the system receives a worklist with the matched RIS procedure name, double click the list, organ programs linked to the matched local procedures directly display on Examination mode for acquisition.

### Procedure for RIS mapping

- ✓ The system is configured for RIS.
- 1 The system receives a DICOM worklist with unknown procedure name.
- 2 Double-click the worklist.

The RIS mapping dialog appears.

- 3 Select category **Adult**, **Paediatric** or **Stitch**, if necessary.
- 4 Click the body part corresponding to the RIS procedure.
- 5 Double-click one or more local procedures which shall be matched to the RIS procedure.

The organ programs assigned to the procedures are shown in the **OGP List**. You can delete one or more organ programs if necessary.

- 6 Click **Match**.

The following message appears: <"Do you want to continue? If 'Yes' is selected, the mapping result will be saved. If 'No' is selected, the mapping result will not be saved.">

- 7 Confirm the match accuracy by clicking **Yes**.

– or –

Click **No** to reselect the local procedure(s).

The following message appears: <"All studies have been mapped. Click 'Yes' to enter examination. Click 'No' to turn back to pre-register page.">

- 8 If you want to continue performing acquisition, click **Yes** to enter **Examination** mode.

– or –

Click **No** to go back to **Patient** mode.


The system establishes the mapping between the procedure name from RIS and the local procedure name permanently. The RIS mapping for the same RIS procedure name need not be performed again.

## Managing RIS mappings

You can modify/delete any mismatching pairs or import a mapping file and perform batch mapping.

- ✓ **Service mode**
- ◆ **Click RIS mapping.**

- Auto RIS mapping**
- 1 Select **Auto RIS mapping** to activate the function.
  - 2 Select the required **RIS mapping DICOM tag**.
  - 3 Click **Apply**.

The following message appears: <"The DICOM tag used for RIS Mapping will be changed. All previous mapping results will not be applicable.">

- 4 Click **Yes** to save the change or click **No** to cancel the change.

- Deleting RIS mappings**
- 1 Select **Adult**, **Paediatric** or **Stitch**.
  - 2 Select the body region by clicking in the figure.

Sort alphabetically

Ascending sort order

Descending sort order

4 Select the procedure.

5 Click **Delete**.

The mapping between the RIS procedure and local procedure(s) is deleted. The RIS procedure name becomes unknown again.

#### Batch mapping

For matching names of local procedures to names of the RIS, it is possible to import a .csv type file, containing RIS names and perform batch mapping for all RIS procedures.

1 Prepare the .csv type file.

The .csv file must contain the following RIS procedure information:

- Procedure code: the value of the DICOM tag used for mapping
- Procedure name: the name of the RIS procedure
- Procedure category: adult/pediatric


It is the responsibility of the operator to provide an accurate list containing the necessary information of RIS procedures to Siemens Healthineers Service to prepare the .csv file.

2 Store the .csv file, e.g. on a USB memory stick.

3 Plug the USB memory into the USB jack.

Wait a moment.

## 4 Click Batch Mapping.

## 5 Click Open and select a mapping file.

The RIS procedure list will populate the batch mapping dialogue box.

## 6 Select the first one to be matched by clicking the list.

## 7 Select Adult, Paediatric or Stitch.

## 8 Select the body region by clicking in the figure.

## 9 Change the sort order of procedures by clicking an icon, if necessary.

By Order

Custom sort order

Ascending sort order

Descending sort order

Sort alphabetically

## 10 Select one or more procedures.

The assigned organ programs are shown in the OGP List.

To delete an organ program from the procedure:

## 11 Select the organ program and click this icon.

To delete all the organ programs from the procedure:

- 12 Click this icon.

To assign the mapping:

- 13 Click **Mapping**.

A confirmation dialog is displayed.

- 14 Confirm the match accuracy by clicking **OK**.

– or –

Click **No** to reselect the local procedures.

- 15 Repeat above steps for all RIS procedures until the batch mapping dialogue box is empty.

#### Auto mapping

It is possible to perform the RIS mapping offline by the Siemens Healthineers service and import the result.

- 1 Prepare the .xml type file.

The .xml file must contain the following RIS procedure information:

- Procedure code: the value of the DICOM tag used for mapping
- Procedure name: the name of the RIS procedure
- Procedure category: adult/pediatric

It is the responsibility of the operator to provide an accurate list containing the necessary information of RIS procedures to Siemens Healthineers Service to prepare the .xml file.

- 2 Click **Auto mapping** and select a mapping file.

Click **Open** and the auto mapping is started.

### 5.4.5 Statistics

You can perform statistical evaluations about your imaging system.

- ✓ **Service mode**

- 1 Click **Statistics**.

- 

**You can print or export the results of FD Battery, Film Supplies, Body Part, and Dosage in \*csv format to a USB drive.**


- ## Audit trail

You can view and export the audit trail.

✓ You are logged in with a user account which has the **Audit Trail** permission.

✓ **Service mode**

1 Click **Audit Trail**.

2 Select **Centralized** and enter the required information of a FTP server.

Click **Echo** to test the connection.

Click **Save** to send the audit trail.

– or –

Select **Locally** to view, remove or export the audit trial to a CD/DVD.

### 5.4.7 Test images

Several test images are available for test purposes.

✓ **Service mode**

◆ Click **Test Image**.

### Displaying the test image overview

- ◆ Click this icon.

### Paging images

To display the next/previous image:

- ◆ Click one of these icons.

### Sending a test image


- ✓ A sending/archiving target node has been configured by Service.
- ◆ Click this icon.  
The image is sent to the configured destination in the background.

### Printing a test image

- ✓ At least one hardcopy device/printer has been configured by Service.
- 1 Select a **Printer** from the list.
  - 2 Select a **Film Size** from the list.

4. **Senders**


- 1000000

If positive, it is \_\_\_\_\_ to the \_\_\_\_\_ quality \_\_\_\_\_ by the \_\_\_\_\_ state \_\_\_\_\_

**Our Technical Customer Service can then use these images for further analysis**

---

- \_\_\_\_\_

- 1 Press the **Ctrl + Alt + S** key combination.

- 2 Describe the error in the blank of **Error Message**.
- 3 Click this icon to select the duration.

If the duration is not selected, only today's log files will be saved.

- 4 Click **Run SaveLog** to start saving.

Operations are not possible during saving the log file.

#### 5.4.11 Service login

For system configuration, maintenance, and diagnosis, service personal can call up the service software.


##### ✓ Service mode

- 1 Click **PMS (Service Software)**.
- 2 Enter the authentication key and click **Login**.  
The service functions are now available.

#### 5.4.12 Remote assistance

If you need help, you can grant remote access to the desktop of the imaging system.

## ✓ Service mode

1 Click **Remote Assistance**.2 To enable remote desktop access: Click **I agree**.

**Connection ID** and **Password** are shown on the window.

3 Inform the **Connection ID** and **Password** to your service engineer.4 Click **Continue**.

A confirmation dialog appears if your service engineer asks for access to your system.

5 Click **Allow** to grant the access.

### 5.4.13 Remote service access

For remote diagnosis it is necessary to enable remote access to the imaging system.

You can check/change access mode at any time.

The following access modes are possible:

- **Full access:** Service is given full access rights to your system.

Operation on your imaging system is blocked.

- **Limited access:** Service is given limited access, that is, they can call up the error messages and reset passwords.
- **No access** (default setting): Service is not given any access to your system.


You can also access the remote service in the **Service** mode.

(→ Page 173 *Remote assistance*)

#### ✓ Service mode

- 1 Click **Remote access**.

- 2 Activate the required item.

- 3 Click **Apply**.

The status icon indicates the remote access is standby.


The status icon indicates an active remote access.

### 5.4.14 DICOM nodes

- ✓ You are logged in with a user account which has the **Service** permission.
- ✓ **Service** mode
- ◆ Click **DICOM Node**.

### Configuring the worklist node

- 1 Click **MWL node**, if necessary.


- 2 Select the **Modality** to query.
- 3 Enter the node and connection data (**Host name**, **IP address**, **IP port**, **Called AE title**, **Calling AE title**).
- 4 Select the **Worklist destination**.
- 5 Click **Enable**.
- 6 Click **Echo** to test the connection.
- 7 Select the **Worklist fields** and the **Matched fields** for study uniqueness.
- 8 Click **Save**.

### Configuring the MPPS node

- 1 Click **MPPS node**.

- 2 Enter the node and connection data (**Host name**, **IP address**, **IP port**, **Called AE title**, **Calling AE title**).
- 3 Select **Enable**.  
If you do not select **Enable**, MPPS configuration is deactivated.
- 4 Click **Echo** to test the connection.
- 5 Click **Save**.

### Configuring sending/archiving nodes

- 1 Click **Archive nodes**.

To add a sending/archive node:

- 2 Click **New**.
- 3 Enter the node and connection data (Host name, IP address, IP port, Called AE title, Calling AE title).
- 4 Click **Archive node** or **Storage commitment node** as required.
- 5 Click **Enable**.
- 6 Click **Echo** to test the connection.

For storage commitment:


- 7 Click **Enable storage commitment** and select the corresponding AE Title.
- 8 Click **Save**.

To delete a sending/archive node:

- 9 Select a node from the **Node name** list and click **Delete**.

## Configuring hardcopy nodes

- 1 Click **Printer nodes**.

To add a printer node:


- 2 Click **New**.
- 3 Enter the node and connection data (Host name, IP address, IP port, Called AE title, Calling AE title).
- 4 Click **Enable**.
- 5 Click **Echo** to test the connection.
- 6 Select the **Default print size**.
- 7 Click **Default printer**, if appropriate.
- 8 Click **Save**.

To delete a printer node:

- 9 Select a node from the **Node name list** and click **Delete**.

## Configuring local DICOM setting

- 1 Click **Local setting**.

- 2 Enter the connection data ( IP address, Subnet mask, Default gateway).
- 3 Click Apply.

#### 5.4.15 FD manager

For Siemens Healthineers service use only

#### 5.4.16 Sensitivity calibration

The detector should be thermalized before calibration.

Make sure that there are no additional objects in the beam path, e.g. a filter.

✓ You are logged in with a user account which has the **Sensitivity** permission.

✓ **Service mode**

1 Click **Sensitivity**.

The following message appears: <"Please insert 0.6 mm copper filter. Please insert default grid for intended plane. Set SID to grid focal distance.">

2 Insert the 0.6 mm Cu filter.

3 Insert the default grid.

4 Set the tube to SID 115 cm.

5 Center the detector in beam field and collimate the beam field to fully cover the active detector area.

6 Click **OK** to start automatic X-ray release.

The Examination task card is displayed.

- 7 Ensure that nobody is in the examination room and release exposure.

The Postprocessing task card is displayed.

- 8 Click Accept.

The following message appears: <"Sensitivity successfully calibrated! Press 'OK' to return to Patient page.">

- 9 Click OK.

The Patient task card is displayed.

### 5.4.17 QR codes

The Imaging system provides QR codes for accessing teamplay FLeet and PEPconnect.

You can open a service ticket on teamplay FLeet and get access to customer educational materials of the system on PEPconnect.

Information is available for dedicated functions, frequently asked questions about the operation of the system, application handling tips, updated instructions, new knowledge, etc.

Example in Service mode

- ◆ Scan the QR code with your mobile device.

teamplay FLeet code: Direct web address for opening ticket will be converted and opened.

PEPconnect code: Direct web address corresponding to MULTIX Impact C system education materials will be converted and opened.

### 5.4.18 HIPAA

The HIPAA (Health Insurance Portability and Accountability Act) feature provides mechanisms to secure the imaging system, patient data and images from misuse by unauthorized persons.

#### Access control

If no user is logged in, it is not possible to access the imaging system, patient or image data. Only after user authentication it is possible to work with the imaging system.

#### Emergency login

The system allows access in a special mode without authentication in order to enable examination in case of an emergency.

A default Emergency User is implemented in the user management with default values and minimal access rights.

For this special user no password is set.

The access privileges of the Emergency User are restricted.

All patients assigned to the Emergency Group or without a group entry (RIS patients) are visible for the Emergency User. All other patients are not visible.

## Administrator login

The administrator can login, using a separate login name and password granting administrator rights to modify the permissions and privileges of all users.

For user and group management, the standard Windows 10® tools are used.

## Data exchange

Sending and retrieving data via DICOM to or from any DICOM server is only possible if the external system is 'known' to the imaging system.

It is assumed, that the system is in a 'secured network' (behind a firewall). For this reason and because of DICOM supports no possibility of data encryption, the imaging system supports the 'trusted hosts' mechanism and does not support data encryption.

## Audit trail

User and system events of the imaging system are recorded in the audit trail.

It is possible to save the audit trail to a USB or CD-R/DVD when logged-in with a user account which has the **Audit Trail** permission.

## Data integrity

The imaging system checks all the entered and received data for integrity. All mandatory data must be available before the imaging system performs any related function (e.g. new patient).

## Data consistency

The imaging system checks all the entered and received data for consistency. This means: format check (e.g. only integer allowed, or no special characters allowed); range check (min/max values; inserted data should be 'useful').

## 6 Glossary

|                   |                                                                                                                                                                                                                                                  |
|-------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 1-point technique | 1-point technique (with Automatic Exposure Control, AEC): An acquisition is taken with operator specified kV value.                                                                                                                              |
| 2-point technique | An acquisition is taken with operator specified kV and mAs values.                                                                                                                                                                               |
| 3-point technique | An acquisition is taken with operator specified kV, mAs and ms values.                                                                                                                                                                           |
| Clinical EXI      | The <b>Clinical EXI</b> analyzes the mean gray value of the diagnostic relevant part in the image via histogram, without the peaks for direct radiation and collimation.                                                                         |
| CR                | Computer Radiography                                                                                                                                                                                                                             |
| DAP               | Dose Area Product                                                                                                                                                                                                                                |
| Deviation Index   | The <b>Deviation Index (DI)</b> detects the tolerance of the exposure in correlation of the <b>Clinical EXI</b> value. The Deviation Index tells something about under- and overexposure of the diagnostic area.                                 |
| DX                | Digital X-ray                                                                                                                                                                                                                                    |
| Image             | DICOM term: An image can be, e.g. a X-ray acquisition or a fluoroscopic image.                                                                                                                                                                   |
| OGP               | Organ program, set of parameters for an acquisition.                                                                                                                                                                                             |
| Patient           | DICOM term: All examination results, i.e. studies, are uniquely assigned to a patient. Unique assignment is usually achieved with the patient's personal data, i.e. last name, first name, date of birth, sex and patient identification number. |
| Procedure step    | DICOM term: Some radiological information systems (RIS) can divide a study into several procedure steps.                                                                                                                                         |
| RAD               | Digital Radiography (RAD mode)                                                                                                                                                                                                                   |
| Series            | DICOM term: A series consists of one or more images. For example, several images acquired in a series are used to display a dynamic process.                                                                                                     |
| SID               | Source to image distance                                                                                                                                                                                                                         |
| SOD               | Source-Object Distance                                                                                                                                                                                                                           |
| STD               | Source-Table Distance                                                                                                                                                                                                                            |
| Study             | DICOM term: A study can consist of single images and image series.                                                                                                                                                                               |
| TOD               | Table-Object Distance                                                                                                                                                                                                                            |

**1, 2, 3 ...**

- 1-point technique 50
- 1:1 pixel size 82
- 2-point technique 50
- 2D camera
  - adjusting ROI for Ortho 68
- 3-point technique 50

**A**

- Accept 73, 87
  - images 58
- Access control 181
  - OGP database 140
- Access mode
  - OGP database 140
- Acq type
  - OGP parameter 150, 152, 155
- acquisition
  - Ortho 70
- Acquisition
  - releasing 55
  - repeating 60
- Acquisition guidance 47
- Acquisition mode 11, 49
- Acquisition parameters
  - checking/changing 48
  - checking/changing on the TUI 54
- Addenda
  - to the Operator Manual 9
- Address
  - configuration 119
- Administration
  - imaging system 118
- Administrator login 182
- AEC field
  - OGP parameter 150
- Angle
  - deleting 97
  - measuring 96
- Application 11
- Archiving
  - configuration 177
  - images 102
  - studies/images 114
- Archiving Information list 26
- Attaching the detector 37
- Audit trail 169, 181, 182
- Auto delete
  - configuration 131
- Auto RIS mapping 165
- Autofilming
  - configuration 130
- Autosending
  - configuration 129

**B**

- Backup
  - OGP database 141, 156
- Batch mapping
  - RIS mapping 166
- Battery failure
  - image recovery 60
- Blend
  - Composing transition mode 75
- BMP
  - exporting 112
- Bucky wall stand
  - Ortho 64

**C**

- Calibration 97
  - by a known distance 100
  - by geometry 98
  - detector 159
- camera
  - adjusting ROI for Ortho 68
- CD
  - exporting 112
  - importing 110
- Changes
  - saving 73, 87
  - undoing 87
- Clearing
  - study 158
- Clinical EXI 56, 57, 58
- Closing
  - a study 61
- Collimation 43
- Collimator
  - OGP parameter 150
- Comments
  - adding 78
  - configuration 123
  - deleting 79
- Composing
  - adjusting manually 75
  - image display 75
  - manual 74
- Configuration 118
  - OGP database 155
- Contrast 92
- Copper filter 53
- Correcting
  - study data 34
- CR 121
  - importing 110
- Cropping
  - an image 81
- Current reduction 52
- Customer data
  - OGP database 142
- Customer mode
  - OGP Editor 140

**Cut**

- Composing transition mode 75

- Cyber security
  - changing 137

**D**

- Daily control study 39
- Data consistency 181, 182
- Data integrity 181, 182
- Date format
  - configuration 121
- DCM
  - exporting 112
  - importing 110
- Deleting
  - studies 35
- Demo image
  - configuration 124
- Department
  - configuration 119
- Detaching the detector 38
- Detail 92
- detailed image 75
- Detector
  - calibration needed 44
  - checking/changing 39
  - firmware upgrade 39
  - management 36
  - status 21
  - temperature 44
- Detector calibration 159
- Detector overview 39
- Deviation Index 56, 58
- DI 58
- DICOM node
  - configuration 175
- DICOM properties
  - configuration 128
- Digital Radiography (RAD) 11
- Disk management
  - configuration 131
- Displaying
  - exam protocol 102
- Distance
  - deleting 97
  - for calibration 100
  - measuring 96, 100
- Domain
  - login 14
- Dominant 52
- Dose 51
  - OGP parameter 150
- Dose setting 57
- DVD
  - exporting 112
  - importing 110

DX 121  
importing 110

## E

ECM  
exporting 112

Edge enhancement 88, 92

EI 57

Emergency login 15, 181

Emergency registration 31

Encryption 13

Evaluation 94

Exam protocol  
displaying 102

Examination  
configuration 124

Examined Patients list 26

Exp correction  
OGP parameter 150

Exp mode  
OGP parameter 150

Export  
OGP database 141, 156

Exporting  
studies/images 112

exposure  
Ortho 70

Exposure  
releasing 55  
repeating 60

Exposure correction 51

Exposure index 56

Exposure parameters  
checking/changing 48  
checking/changing on the TUI 54

Exposure technique 49, 50

Extended exposure 50, 52

## F

F1 18

FD 11  
calibration 159

FD calibration  
performing 160

FD management 36

FD manager 180

FD reconnect 163

Federal law 6

Film layout  
arranging 104

Film preview 104

Film sheet  
arranging 104  
image size 106  
image view 108

Film size  
OGP Editor 155

Filming  
configuration 178

Filter type  
OGP parameter 150

Fit to segment  
film sheet 106

Fit to window 82  
film sheet 106

Flat detector  
image recovery 60

Flat Detector (FD) 11

Flavor  
OGP parameter 152

Flip  
OGP parameter 153

Flipping  
an image 84

Focus 53  
OGP parameter 150

Free texts  
adding 78  
configuration 123  
deleting 79

## G

General settings  
configuration 121

General status 20

Geometrical effects 98

Gray scale  
inverting 81  
inverting on the film sheet 108

Grid 42  
OGP parameter 150

Guide 47

## H

Hardcopy 103  
configuration 178

Heat unit test  
configuration 124

Help 18

Help system 18

HIPAA 181

HIS 26  
searching 28

Hospital Information System (HIS) 26

## I

Illustrations 10

image  
view position and laterality 86

Image  
accepting or rejecting 58  
adding a laterality label 72, 78  
adding text 78

automatic deletion 131

calibrating 97

calibrating by a known distance 100

calibrating by geometry 98

clearing 158

cropping 81

demo image 124

displaying 109

evaluation 94

exporting 112

flavor 88

flipping 84

hardcopy 103

importing 110

inverting the gray scale 81

marking 101

paging 95, 171

panning 83

postprocessing 76

printing 108

processing 76

providing for service 172

Rad 92

Recycle Bin 116

retrieving 116

rotating 84

sending/archiving 102, 114

size on the film sheet 106

storing 101

test 170

transfer status 24

viewing 95

windowing 80

zooming 82

Image acquisition 41

Image display 76  
configuration 125

Image failures  
safety notes 7

Image flip  
OGP parameter 153

Image label  
OGP parameter 153

image laterality 86  
setting 85

Image laterality  
OGP parameter 153

Image processing  
configuration 122

image processing parameter template  
managing 148

Image processing parameters 152

Image receptor size  
configuration 122

Image recovery  
from flat detector 60

Image rotation  
OGP parameter 153

Image text  
configuration 125  
on the film sheet 108

image transition  
   adjusting manually 75  
   display 75

Imaging system 11  
   administration 118  
   operation 17

Import  
   of RIS names 166  
   OGP database 141, 156

Importing  
   studies/images 110

Information  
   about this document 9

Input fields  
   configuration 127

Institution  
   configuration 119

Introduction  
   documents 9

Inverting  
   gray scale 81  
   gray scale on the film sheet 108

**J**

Job status 24

JPG  
   exporting 112

**K**

kV 51  
   OGP parameter 150

**L**

L label  
   adding 72, 78  
   deleting 79

Label  
   adding 72, 78  
   deleting 79  
   language 125  
   OGP parameter 153

Language  
   of laterality label 125  
   of OGP Editor 155  
   user interface 9

laterality 86  
   setting 85

Laterality  
   OGP parameter 153

Laterality label  
   adding 72, 78  
   deleting 79  
   language 125

Latitude 92

Layout 94  
   hardcopy 103

left/right assignment 86

LiveNet 181

Local  
   login 14  
   setting 179

log file  
   providing for service 172

Logging in 14  
   in an emergency 15

Logging off 16

Login 13  
   administrator 182  
   emergency 181  
   service 173

Logoff 13

Look-up table 92

LUT 92

**M**

Magnifying glass 82

Maintenance 173

Manual registration 31

Marking images 101

mAs 51  
   OGP parameter 150

Max. EXI 58

Measurement  
   deleting 97

Measurements 94

Measuring  
   angles 96  
   distances 96, 100

Measuring field 52

Messages 23

Min. EXI 58

Mirroring  
   an image 84

Modality  
   configuration 121, 128

Modality Performed Procedure Step (MPPS) 63

MPPS  
   configuration 128, 176

MPPS close 63

ms 51  
   OGP parameter 150

Multi study mode 26  
   configuration 121

MWL  
   configuration 175

**N**

Name  
   display format 127

Names  
   of patients 10

Noise 92

**O**

Object height  
   calibration 98

OGP  
   checking/changing 45  
   resetting 54  
   RIS mapping 163  
   saving 54

OGP database 141  
   backup and restore 156  
   version, backup 155

OGP Editor 140, 143  
   change password 155  
   changing the password 155  
   closing 157  
   opening 142

Online help 18

Online Help  
   accessing 19

Operation  
   imaging system 17

Operator Manual 9  
   accessing 18  
   scope of applicability 6

Options 9

organ program  
   for Ortho 65

Organ program  
   assigning to procedures 146  
   checking/changing 45  
   configuration 124  
   editing 146  
   editor 140  
   image processing parameters 152  
   miscellaneous parameters 153  
   OGP database 142  
   parameters 149  
   resetting 54  
   RIS mapping 163  
   saving 54  
   selecting 32  
   X-ray parameters 150

Orientation  
   OGP parameter 153

Ortho 64  
   adjusting ROI with camera 68  
   adjusting ROI with start and end position 67  
   composing 71  
   Composing 73  
   on patient table 64  
   on wall stand 64  
   organ program 65  
   positioning on the Ortho stand 66  
   positioning patient on patient table 66  
   preparing 65  
   releasing exposures 70  
   setting 132

overview image 75

**P**

- Paging 95, 171
- Panning
  - an image 83
- Parallax phenomena 98
- Parameter set 92
- Parameters
  - organ programs 149
- Password
  - changing 134, 135
  - changing for OGP Editor 155
  - OGP Editor 155
  - resetting 137
- Patient
  - positioning on patient table 66
  - selecting 30
- Patient administration 109
- Patient data
  - correcting 34
- Patient ID
  - configuration 127
- Patient list
  - configuration 121
  - status 22
- Patient lists 26
- Patient mode 26
- Patient name
  - display format 127
- patient position 86
  - setting 85
- Patient registration 26, 31
  - configuration 127
  - correcting 34
  - emergency 31
  - input fields 127
  - manual 31
  - procedure information 32
- Performed protocol codes 63
- Permissions
  - assigning 136
  - setting up 135
- PMS 173
- Positioning
  - patient on patient table 66
- Postprocessing 76
- Power off 16
- Pre-registered Patients list 26
  - configuration 175
  - updating 28
- Predefined texts
  - adding 78
  - configuration 123
  - deleting 79
- Preview image 45
- Print 103
- Printing
  - configuration 178
  - images 108

- study list 26
- test image 171

- privileges 134

**Procedure**

- data 144
- editing 143
- editor 140
- OGP database 142
- organ programs 146
- RIS mapping 163

- Procedure information 32

- Procedure order
  - configuration 124

- Protection
  - of studies 36

**Q**

- QR codes 181

- Quality control study 39

**R****R label**

- adding 72, 78
- deleting 79

**R/L label**

- language 125

- Rad 92

- RAD 11

- RAD images 11

- Radiological Information System (RIS) 26

- reconstruction algorithm 74

- Recycle Bin 116

- configuration 121

- Recycle Bin list 26

**Registration**

- correcting 34
- emergency 31
- manual 31
- procedure information 32

- Reject 73, 87

- images 58

- Reject reasons

- configuration 124
- entering 58

- Remote access 174

- Remote assistance 173

- Remote service 174

- Remote software update 8

**Reset**

- changes 87
- of OGP settings 54

**Restore**

- OGP database 156

**Retrieving**

- studies/images 116

- RIS 26

- searching 28

- RIS mapping 163

- automatic 165
- batch mapping 166
- deleting 165
- managing 165
- procedure 163

- RIS names

- import 166

**ROI**

- for Ortho with camera 68
- for Ortho with start and end position 67

- roles 134

**Roles**

- adding/changing 136
- assigning 135
- setting up 135

**Rotating**

- an image 84
- on the film sheet 108

**Rotation**

- OGP parameter 153

**S**

- Safety notes 6

- general 6
- image failures 7

- Save 73, 87

- Save As 73, 87

**Saving**

- image 101
- of OGP settings 54

- Scope of applicability 6

- Screen layout 17

**Screen saver**

- configuration 121

**Scrolling**

- in the study list 26

**Searching**

- for studies/patients 28

**Security**

- changing 137

**Selecting**

- a study /patient 30

- Send 94

**Sending**

- configuration 177
- images 102
- studies/images 114
- test image 171

- Sensitivity 51, 57

- Sensitivity calibration 180

- Service functions 158

- Service login 173

- Service mode 173

- Service study 172

- Settings 118

- Shutdown 12, 16

- Shutter 81
  - panning 83
- SID 43, 98
- Siemens data
  - OGP database 142
- Siemens Healthineers mode
  - OGP Editor 140
- Single study mode 26
  - configuration 121
- Site information
  - configuration 119
- Smart virtual Ortho 68
- Smart Virtual Ortho
  - SmartOrtho 64
- SmartOrtho 67
- Smoothing 88, 92
- SOD 98
- Software update 8
- Source to image distance (SID) 43
- Starting 12
- Startup 12
- Statistics 168
- Status
  - of transfers 24
  - patient list 22
- Status information 19
- STD 98
- Stitch
  - mode 71
- Storage commitment
  - configuration 128
- Storing
  - Image 101
- Study
  - clearing 158
  - closing 61
  - deleting 35
  - displaying images 109
  - exporting 112
  - importing 110
  - MPPS close 63
  - protecting/unprotecting 36
  - Recycle Bin 116
  - retrieving 116
  - selecting 30
  - sending/archiving 114
- Study administration 109
- Study data
  - correcting 34
- Study list
  - printing 26
- Study lists 26
  - accessing 26
  - single or multi study mode 26
- Switching off 16
- Switching on 12
- System administration 118
- System messages 23
- System Owner Manual 9
- T
- Target EXI 58
  - OGP parameter 153
- Technical EXI 56
- Temperature
  - of detector 44
- Test image
  - printing 171
  - sending 171
- Test images 170
- Text size
  - changing 78
- Texts
  - adding 78
  - configuration 123
  - deleting 79
  - text size 78
- TIF
  - exporting 112
- Time format
  - configuration 121
- TOD 98
  - for Ortho 67
  - OGP parameter 150
- Top alignment 53
- Transfer 94
  - status display 24
- True size
  - film sheet 106
- Tube load 53, 172
  - OGP parameter 150
- Tube voltage 51
- U
- Undo
  - changes 87
- Unit
  - of length 155
- Unmarking images 101
- Updating
  - the Pre-registered Patients list 28
- Usability 6
- USB
  - configuration 121
  - exporting 112
  - importing 110
- User
  - roles and privileges 134
- User account
  - adding/changing 135
- User interface
  - languages 9
- User management 134
- User settings 118
- Users
  - setting up 135
- V
- Value statements 10
- Version 6
  - OGP Editor, OGP database 155
- View
  - Ortho image 71
  - view position 86
  - setting 85
- View position
  - OGP parameter 153
- Viewing 95
- Virtual film sheet 104
- Virus patterns 8
- W
- Ward
  - configuration 119
- WC
  - shift 88, 92
- Windowing 88, 92
  - an image 80
  - principle 79
  - reset 87
- Wireless detector
  - management 36
  - status 21
- Work list 26
- Worklist
  - configuration 128, 175
  - updating 28
- Workplace 49
- WW
  - shift 88, 92
- X
- X-ray parameter template
  - managing 147
- X-ray parameters 150
- Z
- Zooming
  - an image 82
  - on the film sheet 106

This document is maintained in electronic form through Siemens Healthineers distribution channels. Printed versions not provided by Siemens Healthineers and/or downloaded copies are considered uncontrolled.

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Siemens Shanghai Medical Equipment Ltd. 278 Zhou Zhu Road, Pu  
Dong New Area, Shanghai, P.R. China

Name  
Department

Cheng Xiao  
SHS DI XP QT MEQ

To the Federal Service for Surveillance  
in Healthcare of Russian Federation

Mobile  
E-mail  
Date

+86 (21) 20606826  
cheng.xiao@siemens-healthineers.com  
May 25, 2023

### Self-Truthfulness Declaration for Product Registration

**Legal Manufacturer:** Siemens Shanghai Medical Equipment Ltd.

**Product:** MULTIX Impact C

Herewith we, Siemens Shanghai Medical Equipment Ltd., (SSME), is an important R&D and manufacturing base of Siemens Healthineers covering Computed Tomography (CT), X-Ray products, Mobile C-arms, Ultrasound systems, Digital & Automation as well as Customer Services.

And We, Siemens Shanghai Medical Equipment Ltd., 278 Zhou Zhu Road, Bldg. 2, 3, 4, 5, 7, Pudong New Area, Shanghai, P.R. of China is the legal manufacturer of the following medical device registered on the territory of Russian Federation:

#### X-ray Radiography System MULTIX Impact C with accessories

Registration license No. RZN 2022/18403 dd. 29.09.2022

would like to confirm that the following documents is true and correct copy of the original document:

| Documents                         | Description                                                                                |
|-----------------------------------|--------------------------------------------------------------------------------------------|
| User Manual Supplement for Russia | This document is a supplement manual for user according to Russian regulatory requirements |

Sincerely yours,  
Siemens Shanghai Medical Equipment Ltd.

Wu Jin Yu

Director of SSME XP QT

RECEIVED  
COUNCIL

# User manual supplement for Russia

---

Siemens Shanghai Medical Equipment Ltd.  
278 Zhou Zhu Road, Bldg. 2, 3, 4, 5, 7, Pudong New Area  
Shanghai, P.R. China

## X-ray Radiography System MULTIX Impact C with accessories

## User manual supplement for Russia

### 1. The name of the Medical Device

X-ray Radiography System MULTIX Impact C with accessories

### 2. Information about Manufacturer of the Medical Device:

#### 2.1. Designated name of the manufacturer

##### Legal manufacturer:

Siemens Shanghai Medical Equipment Ltd.  
278 Zhou Zhu Road, Bldg. 2, 3, 4, 5, 7, Pudong New Area  
Shanghai, P.R. China  
Tel.: +86 21 3889 5000  
Email: [www.siemens-healthineers.com](http://www.siemens-healthineers.com)

#### 2.2. Addresses of the additional manufacturing sites;

##### Manufacturing facility:

1. Siemens Shanghai Medical Equipment Ltd. XP Facility  
278 Zhou Zhu Road, Shanghai 201318, P.R.China
2. Siemens Healthcare GmbH, X-ray products  
Siemensstr.1, 91301 Forchheim, Germany

#### 2.3. The name and address of the authorized representative in Russia.

##### Authorized representative in Russia:

Siemens Healthcare LLC  
Dubinskaja str. 96, 115093 Moscow, Russia  
Tel: (495) 737 12 52

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#### 3.2. Intended purpose;

*Intended Purpose Statement under the European Medical Device Regulation 2017/745:*

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The climatic conditions are specified below:

| Examination room    |                            |
|---------------------|----------------------------|
| Temperature range   | +10°C to +35°C             |
| Relative humidity   | 20% to 75%, non-condensing |
| Barometric pressure | 700 hPa to 1060 hPa        |
| Height              | ≤ 3000 m above sea level   |
| Imaging system      |                            |
| Temperature range   | 0°C to +35°C               |
| Relative humidity   | 20% to 75%, non-condensing |

## 4. Classification of the Medical Device;

### 4.1. Type of the MD in accordance with the nomenclature for classification;

Kind of product: Active medical device for diagnosis

Classification rules according to Medical Device Regulation (EU) 2017/745, Annex VIII

Classification according to Annex VIII: IIb

### 4.2. Class of the potential risk of the medical device application;

SW safety classification acc. to EN 62304: based on the current product risk management results, the SW safety class – is "B"

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*Attention: To avoid the risk of electrical shock, a protective conductor must be implemented.*

### 4.4. Applied parts;

- Table
- Compression device
- Wall stand
- Portable detector

### 4.5. Protection against ingress of water;

IPx8: Foot switches

IP43, IP67: MAX wi-D

IPx3: Core XL

IPx0: Rest of the system

#### 4.6. Explosion protection (AP, APG after IEC 60601-1-3:2013; IEC 60601-1:2012);

The system is not intended for operation in explosion-endangered areas (for example highly flammable mixtures of anaesthesia gases with air or oxygen or nitrous oxide).

#### 4.7. Operating mode;

Continuous operating with intermitted loading.

#### 4.8. Degree of contamination;

Class 2, according to IEC 60601-1-3:2013; IEC 60601-1:2012.

#### 4.9. Laser classification;

For Collimator AL04:

| Laser localizer |        |
|-----------------|--------|
| Class           | 2      |
| Wave length     | 650 nm |
| Power           | <1mW   |

For Collimator RFU:

| Laser localizer |        |
|-----------------|--------|
| Class           | 1      |
| Wave length     | 650 nm |
| Power           | <1mW   |

The lasers are the same for the two collimators. The laser classes are different as the power of two collimators are different.

### 5. Technical Specifications of the Medical Device (the main technical parameters of the Medical Device; the MD's scheme (design), characteristics, properties, features):

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- On the patient table, X-rays can be taken on the recumbent or sitting patient from head to foot. Exposures in the region of the skull, of the spine, of thorax, lungs and abdomen as well as of the extremities are possible.

- On the wall stand, X-rays can be taken on the standing or sitting patient (chest, and so on).


In addition, exposures onto the wireless detectors in the form of free exposures are possible.

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



- MULTIX Impact C with ceiling-mounted tube support (ACSS & Tracking):


- Bucky wall stand with Core static
- Bucky wall stand with Core static, and optional Core XL
- Bucky wall stand with Core static, and optional MAX wi-D
- Bucky wall stand with Core XL
- Bucky wall stand with MAX wi-D
- Bucky wall stand with Core static and patient table with Core XL
- Bucky wall stand with Core static and patient table with MAX wi-D
- Bucky wall stand and patient table to be used with one Core XL
- Bucky wall stand and patient table to be used with one MAX wi-D

Please see summary below:

|   | System solutions                                                                                                                                                           | Patient table | Detector for patient table | BWS                       | Detector for BWS |
|---|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---------------|----------------------------|---------------------------|------------------|
| 1 |  <p><i>Ceiling-mounted solution with Bucky wall stand equipped with Core static</i></p> | -             | -                          | BWS (integrated detector) | Core static      |

|   |                                                                                                                                                                                                       |   |   |                                                    |                         |
|---|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|---|---|----------------------------------------------------|-------------------------|
| 2 |  <p><i>Ceiling-mounted solution with Bucky wall stand equipped with Core static and Core XL for free exams</i></p>   | - | - | BWS<br>(integrated detector and wireless detector) | Core static<br>Core XL  |
| 3 |  <p><i>Ceiling-mounted solution with Bucky wall stand equipped with Core static and MAX wi-D for free exams</i></p> | - | - | BWS<br>(integrated detector and wireless detector) | Core static<br>MAX wi-D |
| 4 |  <p><i>Ceiling-mounted solution with Bucky wall stand equipped with Core XL</i></p>                                | - | - | BWS<br>(wireless detector)                         | Core XL                 |
| 5 |  <p><i>Ceiling-mounted solution with Bucky wall stand equipped with MAX wi-D</i></p>                               | - | - | BWS<br>(wireless detector)                         | MAX wi-D                |

Ceiling-mounted solution with Bucky wall stand equipped with Core static as well as height adjustable table equipped with Core XL


Available  
(wireless  
detector)

Core XL

BWS  
(integrated  
detector)

Core static

7

Ceiling-mounted solution with Bucky wall stand equipped with Core static as well as height adjustable table equipped with MAX wi-D


Available  
(wireless  
detector)

MAX wi-  
D

BWS  
(integrated  
detector)

Core static

8

Ceiling-mounted solution with Bucky wall stand and height adjustable table.  
Core XL can be used either in the table or Bucky wall stand


Available  
(wireless  
detector)

Core XL\*

BWS  
(wireless  
detector)

Core XL

9

Ceiling-mounted solution with Bucky wall stand and height adjustable table.

Available  
(wireless  
detector)

MAX wi-  
D\*

BWS  
(wireless  
detector)

MAX wi-D

|                                                              |  |  |  |  |
|--------------------------------------------------------------|--|--|--|--|
| MAX wi-D can be used either in the table or Bucky wall stand |  |  |  |  |
|--------------------------------------------------------------|--|--|--|--|

*\*BWS and Patient table can be used with one wireless detector*

The list of components covers all the possible configurations. Please kindly find a basic blocks summary below:

| Basic block components                            | Basic block 1            | Basic block 2            |
|---------------------------------------------------|--------------------------|--------------------------|
| Tube stand 3D: standard or motorized              | Available                | Available                |
| Bucky wall stand with Core static                 | Available (if necessary) | Available (if necessary) |
| Bucky wall stand for MAX wi-D                     | Available (if necessary) | Available (if necessary) |
| Bucky wall stand for Core XL                      | Available (if necessary) | Available (if necessary) |
| MAX wi-D detector                                 | Available (if necessary) | Available (if necessary) |
| Core XL detector                                  | Available (if necessary) | Available (if necessary) |
| X-ray generator:<br>- 55 kW<br>- 65 kW<br>- 80 kW | Available                | Available                |
| IS workstation w/touch-screen                     | Available (if necessary) | Available (if necessary) |
| IS workstation standard (w/o touch screen)        | Available (if necessary) | Available (if necessary) |
| CAREMAX plus HS Integrated                        | Available                | Available                |
| User documentation                                | Available                | Available                |
| Table for MAX wi-D                                | Not available            | Available (if necessary) |
| Table for Core XL                                 | Not available            | Available (if necessary) |

All the technical specifications and the list of components are described in the chapters below.

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| Installation data                |                                                                   |
|----------------------------------|-------------------------------------------------------------------|
| Installation data                | The entire system is powered via a three phase voltage connection |
| Line voltage connection (system) | 3-phase, 380V, 400V, 440V(50/60Hz)± 10%, 480V(60Hz) ± 10%         |
| Power consumption                | Max. 127kVA(80kW)                                                 |
| Standby mode                     | 1.2 kVA                                                           |

## Internal line resistance

| ULine | 55 kW  | 65 kW  | 80 kW  |
|-------|--------|--------|--------|
| 380 V | 0,15 Ω | 0,15 Ω | 0,10 Ω |
| 400 V | 0,17 Ω | 0,17 Ω | 0,11 Ω |
| 440 V | 0,20 Ω | 0,20 Ω | 0,14 Ω |
| 480 V | 0,24 Ω | 0,24 Ω | 0,16 Ω |

| Weight of the main components |                |
|-------------------------------|----------------|
| Bucky wall stand              | Approx. 240 kg |

|                    |                                                                                      |
|--------------------|--------------------------------------------------------------------------------------|
| X-ray tube support | With 3 m rail: approx. 300 kg<br>With 4 m rail: approx. 320 kg                       |
| Imaging station    | With normal screen: approx. 14,5 kg<br>With touch screen (optional): approx. 15,4 kg |
| Generator cabinet  | approx. 190kg                                                                        |
| Patient table      | approx. 320 kg                                                                       |

### 5.3. The list of components of the medical device;

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|                                                                                                                                                        |                              |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
|--------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| D (при необходимости).                                                                                                                                 |                              | <p>Detector tray can be operated from the left or right side (depends on the Bucky wall stand configuration). Vertical height adjustment is possible from left and right side. The vertical height adjustment can be also performed manually.</p> <p>It is especially suited for acquisitions of skeletal radiography of the standing and seated patient:</p> <ul style="list-style-type: none"> <li>- Orthopedic diagnostics</li> <li>- Thorax and general diagnostics</li> <li>- Trauma and ER diagnostics.</li> </ul>                                                                                                                                                                               |
| Стойка вертикальная для беспроводного детектора Core XL (при необходимости).                                                                           | Bucky wall stand for Core XL | <p>Floor - mounted Bucky wall stand with motorized height - adjustable detector tray for a Core XL 43 cm x 43 cm (17" x 17") wireless detector for digital acquisitions.</p> <p>Detector tray can be operated from the left or right side (depends on the Bucky wall stand configuration). Vertical height adjustment is possible from left and right side. The vertical height - adjustment can also be performed manually.</p> <p>It is especially suited for acquisitions of skeletal radiography of the standing and seated patient:</p> <ul style="list-style-type: none"> <li>- Orthopedic diagnostics</li> <li>- Thorax and general diagnostics</li> <li>- Trauma and ER diagnostics</li> </ul> |
| Стол пациента для детектора MAX wi-D (при необходимости).                                                                                              | Table for MAX wi-D           | Height - adjustable patient positioning table in compact design, for X-ray exposures of the entire body with detector tray for MAX wi-D 35 cm x 43 cm (14" x 17").                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Стол пациента для детектора Core XL (при необходимости).                                                                                               | Table for Core XL            | Height - adjustable patient positioning table in compact design, for X-ray exposures of the entire body with floating tabletop and detector tray for Core XL 43 cm x 43 cm (17" x 17") wireless detector.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Беспроводной плоский детектор MAX wi-D (не более 5 шт., при необходимости).                                                                            | MAX wi-D                     | <p>Light, wireless, mobile detector 35 cm x 43 cm (14" x 17"), with handle for comfortable and safe operation. Thanks to the MAXswap function, the detector can be used with all MAX systems.</p> <p>It can be automatically charged in the system's detector holder.</p>                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Беспроводной плоский детектор Core XL (не более 5 шт., при необходимости).                                                                             | Core XL                      | <p>Light - weight mobile, wireless 43 cm x 43 cm (17" x 17") detector with handle for comfortable and safe handling.</p> <p>It can be charged automatically in the system's detector holder.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <p>Генератор рентгеновский, варианты исполнения:</p> <ul style="list-style-type: none"> <li>- 55 кВт;</li> <li>- 65 кВт;</li> <li>- 80 кВт.</li> </ul> | 80 kW generator              | <p>High-frequency X-ray generator with 80 kW power, with high performance to extend the range of applications.</p> <p>Performance data:</p> <ul style="list-style-type: none"> <li>- Max. power: 80 kW at 100 kV</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |

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|                                                                      |                                   |                                                                                                                                                                                                                                                                                                                                                                                                |
|----------------------------------------------------------------------|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                      |                                   | <p>- Grid focusing for source-image distance (SID) of 140 cm (55")</p> <p>- Universal grid: Working range between F115 cm and F180 cm (45" to 71")</p> <p>- Dimensions without grip (W x H): 460 mm x 460 mm (18.11" x 18.11")</p> <p>Weight: 1.30 kg (2.86 lbs)</p>                                                                                                                           |
| Решетка отсеивающая (13/92, F180) (не более 5 шт.).                  | Transparent grid 13/92 F180       | Highly selective anti-scatter grid for scattered radiation reduction.                                                                                                                                                                                                                                                                                                                          |
| Решетка отсеивающая (13/40, F115) (не более 5 шт.).                  | Transparent grid 13/40 F115       | Highly selective anti-scatter grid for scattered radiation reduction.                                                                                                                                                                                                                                                                                                                          |
| Решетка отсеивающая универсальная (13/40) (не более 5 шт.).          | Transparent grid 13/40, universal | Highly selective anti-scatter grid for scattered radiation reduction.                                                                                                                                                                                                                                                                                                                          |
| Решетка отсеивающая (13/40, F180) (не более 5 шт.).                  | Transparent grid 13/40 F180       | Highly selective anti-scatter grid for scattered radiation reduction.                                                                                                                                                                                                                                                                                                                          |
| Решетка отсеивающая (13/92, F300) (не более 5 шт.).                  | Transparent grid 13/92, F300      | Highly selective anti-scatter grid for scattered radiation reduction.                                                                                                                                                                                                                                                                                                                          |
| Зарядное устройство для батарей детектора MAX wi-D (не более 5 шт.). | Detector battery charger MAX wi-D | <p>This charger can be used to charge the replacement batteries for the MAX wi-D detector.</p> <p>Space for 3 batteries, with LED indicator for charge status. The charger connects to a wall socket using a power cord.</p> <p>This item includes the following components:</p> <ul style="list-style-type: none"> <li>- 1x battery charger</li> <li>- 1x detector battery</li> </ul>         |
| Зарядное устройство для батарей детектора Core XL (не более 5 шт.).  | Detector battery charger Core XL  | <p>This charger can be used to charge the replacement batteries for the Core XL wireless detector.</p> <p>Space for one battery, with LED indicator for charge status. The charger connects to a wall socket using a power cord.</p> <p>This item includes the following components:</p> <ul style="list-style-type: none"> <li>- 1x battery charger</li> <li>- 1x detector battery</li> </ul> |
| Точка доступа беспроводного                                          | WLAN for wireless detectors       | Integrated WLAN access point for operation of the MAX wi-D and CORE Core XL detector.                                                                                                                                                                                                                                                                                                          |

|                                                                              |                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------------------------------------------------------------|------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| соединения WLAN для детекторов.                                              |                                    |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
| Камера для позиционирования пациента (не более 5 шт.).                       | myExam 2D camera                   | <p>Camera attached to the collimator streams live images from the Bucky wall stand and/or table to the imaging station. By ensuring proper positioning this may reduce the number of retakes necessary and thereby lessen the dose per requested image.</p> <p>Only available as option and can only be configured with motorized collimator.</p>                                                                                                                                                                                |
| Программный режим изменения мощности генератора.                             | 65 kW generator upgrade            | <p>Upgrade high-frequency X-ray generator to 65 kW power, for increasing performance to extend the range of applications.</p> <p><u>Performance data:</u></p> <p>- Max. power: 65 kW at 100 kV</p> <p>Max. tube current: 800 mA</p>                                                                                                                                                                                                                                                                                              |
| Программный режим автоматического позиционирования системы SmartPositioning. | SmartPositioning                   | <p>SmartPositioning accelerates workflows meaning the system can automatically adjust the height of the tube, table, and Bucky wall stand and the desired angle by tilting the tube automatically. With the push of one button it moves MULTIX Impact C into a pre-defined and OGP-linked positions only necessary to adjust the x- and y-axis of the ceiling-mounted tube manually.</p> <p>This reduces the overall physical burden of system operation and brings more ease into preparing and setting up the examination.</p> |
| Программный режим автоматического отслеживания детектора для стола и стойки. | Table Bucky Auto tracking          | <p>Enable tube &amp; table detector horizontal and oblique tracking. It offers oblique and horizontal tube tracking with table detector for greater accuracy. This function can enhance system automation for faster examinations and less physical burden on staff.</p>                                                                                                                                                                                                                                                         |
| Подставка для станции формирования изображений.                              | Stand f. IS workstation w/touch s. | <p>Desktop stand, as option for image station with optional integrated 23.8" color touch screen.</p>                                                                                                                                                                                                                                                                                                                                                                                                                             |
| Педальный переключатель.                                                     | Foot switch (3 pedals)             | <p>3-pedal foot switch for adjusting the height of the patient table and releasing the floating tabletop. Technical details:</p> <p>- Cable length: 4.5 m (177")</p> <p>- Weight: 2.2 kg (4.8 lbs)</p>                                                                                                                                                                                                                                                                                                                           |
| Ножной переключатель двухсторонний.                                          | Upgr Foot Kick Switch Front & Rear | <p>For height adjustment of the patient positioning table and switching of the floating tabletop. Height-adjustment, release, and locking of the floating tabletop is done through a foot kick switch. The foot kick rails are located in the foot area both</p>                                                                                                                                                                                                                                                                 |

|                                                               |                                  |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|---------------------------------------------------------------|----------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                               |                                  | at the front side and the rear side of the patient positioning table and can be programmed individually at the time of installation. This prevents accidental operation by patients or accompanying persons.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Беспроводной пульт дистанционного управления.                 | Wireless remote control          | <p>Touch sensitive wireless remote control enables convenient control over Bucky wall stand vertical movements and collimation field size.</p> <p>Only available as option and can be only configured with motorized collimator.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Ручной переключатель экспозиции.                              | Exposure hand-switch             | <p>Handswitch for radiation release, with cable connection to the generator ON/OFF console.</p> <p>Only available as option.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Программный режим SmartOrtho.                                 | SmartOrtho                       | <p>SmartOrtho is an automated tilting technique for long-leg and full-spine imaging. Up to 4 single images can be acquired to cover the selected region of interest with the patient in standing or lying position. The images are automatically composed into a single image on the workstation.</p> <ul style="list-style-type: none"> <li>- Acquisitions made at the Bucky wall stand, should use a source-image distance (SID) of 180 cm or 300 cm.</li> <li>- Acquisitions made at the table, should use a source-image distance (SID) of 130 cm.</li> <li>- The use of an 80 kW generator is recommended for acquisitions of the lateral spine.</li> <li>- SmartMove /SmartPositioning is recommended for fast and automatic system movement to the OGP- linked and predefined position.</li> </ul> <p>SmartOrtho License Key</p> |
| Программный режим виртуальной коллимации Smart Virtual Ortho. | Smart Virtual Ortho              | <p>The myExam 2D camera live image enables accurate collimation and helps users to manually adjust ortho planning area at Bucky wall stand and table on the workstation.</p> <p>Smart Virtual Ortho allows you to accurately set the start and end point of the imaging. Patient movements can be detected quickly, avoiding unnecessary retakes, and easing the burden on orthopedic patients. Via Auto Stitching, up to four single images are automatically acquired and processed into a single image, leaving more time for patient care. With Smart Virtual Ortho you can easily judge whether you can perform an exam with fewer images - which can help reduce dose for each individual patient.</p>                                                                                                                            |
| Программный режим                                             | Auto Full-Spine&Long-Leg Collim. | Automatic detection of the full-spine or long-leg and provide ortho planning area for standing                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |

|                                                                    |                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                           |
|--------------------------------------------------------------------|-------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| автоматической коллимации.                                         |                                     | <p>examinations at the Bucky wall stand via the myExam 2D Camera.</p> <p>An AI-based algorithm detects the full spine or long leg and collimates accordingly, leaving your staff to fully focus on their patients.</p> <p>Auto long-leg and Auto full-spine Collimation have built-in Artificial Intelligence that automatically detects the respective region of interest and sets the appropriate collimation automatically. This helps you deliver long-leg and full-spine images with ease. You can focus on your patients and their well-being - while the system takes care of everything else.</p> |
| Лицензия временная на программный режим автоматической коллимации. | Auto Full-Spine&Long-Leg Coll TRIAL | <p>90 days trial of the Auto full-spine &amp; long-leg Collimation.</p> <p>Automatic detection of the full-spine or long-leg and provide ortho planning area for standing examinations at the Bucky wall stand via the myExam 2D Camera.</p> <p>An AI-based algorithm detects the full spine or long leg and collimates accordingly, leaving your staff to fully focus on their patients.</p>                                                                                                                                                                                                             |
| Универсальная стойка Ortho (не более 3 шт.).                       | Multipurpose/Ortho stand            | <p>The multipurpose stand stabilizes the patient during upright examination, such as long -leg and full-spine imaging to prevent motion artifacts.</p> <p>It allows the safe movement of the Bucky tray during the image acquisition procedure. The stand is light weight, easy to maneuver even through doors and provides height adjustable handgrips for the patient.</p>                                                                                                                                                                                                                              |
| Линейка с держателем для стойки Ortho.                             | Multipurpose stand ruler and holder | Ruler for measurements on acquired images (e.g. on composed long leg images). The ruler holder slides onto a rail at the column of the stand and can be adjusted according to region of interest.                                                                                                                                                                                                                                                                                                                                                                                                         |
| Дисковод для записи CD/DVD.                                        | DVD-W ROW                           | CD/DVD Read & Write drive for digital image storage on CD-R/DVD for offline data exchange in DICOM format.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| Программный режим CAREMAX Virtual.                                 | CAREMAX Virtual                     | <p>CAREMAX Virtual is a software solution for acquisition of the dose-area product (DAP). Resolution of 0.01 <math>\mu\text{Gym}^2</math>.</p> <p>Comparing CAREMAX option, CAREMAX Virtual doesn't include any hardware (e.g. electronics unit). The software solution provided by CAREMAX Virtual should be activated via a license.</p> <p>With CAREMAX Virtual the calculated dose-area product is displayed on the flat-screen display.</p>                                                                                                                                                          |

|                                                              |                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
|--------------------------------------------------------------|-------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| осветка<br>диодная.                                          | Impact C Illuminate                 | <p>LED lights on top of the ceiling stand and the Bucky wall stand indicate the status and progress of the imaging - so patients know exactly when they can stop holding their breath, for instance.</p> <p>Only available as option.<br/>3 colors indicate the different system status:</p> <ul style="list-style-type: none"> <li>- LED green: System is ready for exposure</li> <li>LED yellow: Exposures are being performed</li> <li>LED blue: System is not ready for exposure</li> </ul>                                                                                                                                                                                                                                     |
| Устройство для дистанционного управления на основе планшета. | Remote Interface (HW & license key) | Tablet-based workflow for the system. Includes multiple functions, for example, to check the list of registered patients when picking patient up in the waiting room. You can also identify patients by scanning their QR or barcodes before you take them to the examination room and edit the list of organ programs (OGPs) on the go.                                                                                                                                                                                                                                                                                                                                                                                            |
| Удлинение потолочных рельсов.                                | Track Extension                     | Extension of the longitudinal travel range of the ceiling-mounted support by a maximum of 2.12 meters.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                              |
| Удлинение штатива телескопическое.                           | Telescopic extension                | In high rooms the ceiling-mounted tube support can be height-extended by 200 mm.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                    |
| Источник бесперебойного питания (не более 3 шт.).            | UPS Eaton for Imaging System        | <p>Universal Power Supply for the imaging system to ensure safe and proper shutdown in case of a power interruption.</p> <p>Line-interactive high frequency UPS with:</p> <ul style="list-style-type: none"> <li>- Power of 600 W</li> <li>- Size: 230 mm x 150 mm x 345 mm</li> <li>- Weight: 10.4 kg</li> </ul>                                                                                                                                                                                                                                                                                                                                                                                                                   |
| Держатель мобильный.                                         | Detector holder, mobile             | <p>Mobile detector holder for safe placement of the mobile detector, CR and film cassettes from 24 cm x 30 cm to 43 cm x 43 cm (9.5" x 12" to 17" x 17"), e.g. for cross-table or other free examinations.</p> <p>The detector holder can be used with the mobile detector, CR and film cassettes / grids from 24 cm x 30 cm to 43 cm x 43 cm (9.5" x 12" to 17" x 17").</p> <p>Vertical movements are counterbalanced for simple and quick height adjustments. The holder is suitable for cross-table views. The sides of the holder can be quickly extended or retracted to secure different detector / cassette sizes. The holder can be turned, tilted and orientated to suit any examination position infinitely variable.</p> |

|                                                                      |                                     |                                                                                                                                                                                                                                                                                                                                                                                                                                 |
|----------------------------------------------------------------------|-------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                      |                                     | <p>- Dimensions (D x W x H): 109 cm x 61 cm x 179 cm (43" x 24" x 70.5")</p> <p>Height adjustable from floor level to 129 cm (50.8")</p> <p>Overhangs table or trolley by up to 62 cm (24.4")</p> <p>Telescopic width adjustment: 24.5 cm to 53.5 cm (9.7" x 21")</p> <p>Weight: 55 kg (121 lbs)</p>                                                                                                                            |
| Опорная ручка для вертикальной стойки.                               | Patient handle lateral              | For lateral exposure of thorax on wall stand.                                                                                                                                                                                                                                                                                                                                                                                   |
| Держатель детектора для универсальной стойки Ortho (не более 3 шт.). | Detector holder for multipur. stand | Attach a wireless detector to the multipurpose stand e.g. for pediatric imaging. The detector holder slides onto a rail at the column of the stand and can be operated with one hand. It holds and automatically centers the MAX wireless detectors in portrait and landscape orientation.                                                                                                                                      |
| Дополнительный аккумулятор для детекторов (не более 3 шт.).          | Battery MAX mini and MAX wi-D       | <p>Additional battery for MAX wi-D and MAX mini. The battery can be loaded in a Bucky tray while inside the MAX wi-D or in the battery charger.</p> <p>Technical details:</p> <p>- Dimensions: 24.7 cm x 7.4 cm x 0.7 cm (9.7" x 2.9" x 0.28")</p> <p>Weight: 100 g (0.220 lbs)</p> <p>Battery capacity: 3,35 Ah</p> <p>Max. voltage: 7.4V</p>                                                                                  |
| Дополнительный аккумулятор для детектора Core XL (не более 3 шт.).   | Battery Core XL                     | Additional battery for Core XL wireless detector. The battery can be charged in a Bucky tray (BWS or table) while inserted in the Core XL detector or in the optional battery charger.                                                                                                                                                                                                                                          |
| Держатель детектора (боковой).                                       | Detector holder lateral             | <p>Detector holder tray for direct exposures transverse to the table for lateral beam projection and horizontal tabletop. The detector is positioned vertically. A wall bracket is included in the scope of delivery. Detector with and without attachable grid.</p> <p>Suitable for 35 x 43 MAX wi-D detector and CR cassettes.</p> <p>Technical details:</p> <p>- Dimensions (H x W): 22 cm x 53 cm</p> <p>Weight: 5.3 kg</p> |
| Сканер штрихового кода.                                              | MULTIX Impact barcode reader        | The barcode reader is used for scanning patient data which stored in barcodes.                                                                                                                                                                                                                                                                                                                                                  |

|                                   |                               |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
|-----------------------------------|-------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                   |                               | Usage example in patient mode: the barcode reader can be used for scanning patient ID and helps for efficient searching and matching the pre-registered and/or examined patients' list.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Фиксирующий ремень.               | Compression belt              | <p>Belt compression device used for patient compression during thoracic or abdominal examinations as well as for safe positioning of restless or frail patients. Compression is achieved by means of an easy to clean, radiolucent plastic belt with a ratcheting tension lever.</p> <p>Advantages of compression:</p> <ul style="list-style-type: none"> <li>– Quick and safe securing of patient to the tabletop.</li> <li>– Reduction of patient thickness, i.e. improvement of image quality through reduction of scattered radiation.</li> </ul> <p>Technical details:</p> <ul style="list-style-type: none"> <li>- Weight: 3,15 kg</li> <li>- Dimensions: 16 cm x 36 cm x 12 cm (L x B x H)</li> <li>- Belt width: 23 cm</li> <li>- Belt length (completely unrolled): 144 cm</li> </ul> |
| Держатель кассеты без фиксатора.  | Cassette holder without clamp | <p>For cross-table direct exposures, lateral beam path, and horizontal tabletop.</p> <p>Technical details:</p> <p>The cassette holder supports the cassette in an upright position anywhere on the horizontal tabletop according to the desired cross-table exposure. For IEC, DIN and ANSI-format cassettes from 13 cm x 18 cm (5" x 7") to 35 cm x 43 cm (14" x 17"), vertical and horizontal.</p> <p>Weight: 1.5 kg (3.3 lbs)</p>                                                                                                                                                                                                                                                                                                                                                           |
| Планка-поручень (не более 5 шт.). | Hand grip rail waved          | <p>Detachable handgrip rail (waved) that can be attached to the lateral tabletop accessory rails for safe patient positioning.</p> <p>Standard item for fluoroscopy systems, included in basic unit and supplied with system. Only available as an option for RAD systems. Also optionally available as a spare part.</p> <p>Length: 74 cm</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| Опорная ручка.                    | Hand grip oblique             | <p>Handgrips enable the patient to hold on during the examination.</p> <p>They are mounted directly to the tabletop using clamping screws.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |

|                                                           |                                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
|-----------------------------------------------------------|----------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                           |                                        | Weight: 0.5 kg                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| Матрац пациента.                                          | Patient positioning mattress           | <p>Radiolucent table pad with a heavy-duty, soft, light-colored plastic cover that is easy to clean. The soft cushion allows comfortable patient positioning and repositioning. To prevent the pad from sliding during head-up positions, the straps of the patient table pad can be attached to the grip protection rail at the head end.</p> <p>The soft cushion allows comfortable patient positioning and repositioning.</p> <p>Technical details:</p> <ul style="list-style-type: none"> <li>- Length: 198 cm (78")</li> <li>- Width: 66 cm/25.9" (of which 53.5 cm is padded)</li> <li>- Thickness: 2.5 cm (0.98")</li> <li>- Weight: 2.7 kg (5.9 lbs)</li> </ul> |
| Вращающийся стул пациента для универсальной стойки Ortho. | Patient stool<br>f.Multipurp./3D stand | <p>Lightweight, height adjustable patient stool to easily position patients (e.g. pediatric, geriatric). The attachment of the stool to the stand is mechanically guided to minimize positioning efforts. The sitting height adjustment can be done with the operator's foot so the hands are free to support the patient. The distance between stool and screen can also be adjusted. Because the stool can rotate 360°, there is no need for the patient to stand up and be repositioned during acquisitions while sitting (e.g. shoulder).</p>                                                                                                                       |
| Держатель рулона бумаги.                                  | Paper roll holder                      | <p>The paper roll holder is attached to the end of the table; the paper roll is inserted slightly below the table height. Paper rolls of different diameters and lengths may be used.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                               |
| Стол для панели управления (не более 3 шт.).              | High desk, fixed height                | <p>Computer desk for the image and control station.</p> <p>Technical details:</p> <ul style="list-style-type: none"> <li>- Working height of 100 cm (39.4")</li> </ul> <p>Desktop working area of 80 cm (31") x 100 cm (39.4")</p>                                                                                                                                                                                                                                                                                                                                                                                                                                      |
| Настенный держатель для решеток (не более 3 шт.).         | Wall holder for grid                   | <p>Wall bracket with 2 compartments of different widths to store replaceable grids or cassette drawers (weight: 5.0 kg/11 lbs).</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| Марки для разметки пленки (не более 10 комплектов).       | Film marking R/L                       | <p>Film marks for marking of right and left information on film. The mark is with a magnet holder.</p>                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |

|                                                             |                                 |                                                                                                                                                                                                                                                                                                                                                                                                                   |
|-------------------------------------------------------------|---------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Интеркомное устройство врач-рентгенолог.                    | Intercom System Basic           | Intercom system for communication between examination room and control room.                                                                                                                                                                                                                                                                                                                                      |
| Держатель BABIX (не более 3 шт.).                           | Baby carrier holder BWS         | The baby carrier holder is used to fasten carrier for examination to the Bucky wall stand during pediatric examinations of babies. The baby carrier is not part of the delivery.                                                                                                                                                                                                                                  |
| Держатель детской люльки для стойки Ortho (не более 3 шт.). | Baby carrier holder Ortho Stand | The baby carrier holder is used to fasten carrier for examination to the multipurpose / ortho stand during pediatric examinations of babies. The baby carrier is not part of the delivery.                                                                                                                                                                                                                        |
| Компенсационный фильтр, таз (не более 3 шт.).               | Homogenizing filter pelvis      | For compensation of density variations for pelvis exposures. Compensating filters prevent overexposures and contribute significantly to improved image quality.<br>The compensation filter ensures object-adapted absorption compensation for bones / tissue or air transitions. They are made of aluminum and can be attached to the multileaf collimator.<br><br>Weight: 0.6 kg (1.32 lbs)                      |
| Компенсационный фильтр, череп взрослого (не более 3 шт.).   | Homogenizing filter adult skull | For compensation of density variations for exposures of the head and neck region of adults. Compensation filters prevent overexposure and contribute significantly to improved image quality. The aluminum compensation filter provides object-specific absorption compensation for bone, tissue and air transitions. They can be attached to the multileaf collimator.<br><br>Weight: 0.6 kg (1.32 lbs)          |
| Компенсационный фильтр, нога (не более 3 шт.).              | Homogenizing filter foot        | For compensation of density variations for exposures of the foot region. Compensation filters prevent overexposure and contribute significantly to improved image quality.<br>The aluminum compensation filter provides object-specific absorption compensation for bone, tissue and air transitions. They are made of aluminum and can be attached to the multileaf collimator.<br><br>Weight: 0.6 kg (1.32 lbs) |
| Компенсационный фильтр, череп ребенка (не более 3 шт.).     | Homogenizing filter child skull | For compensation of density variations for exposures of the head and neck region of children. Compensation filters prevent overexposure and contribute significantly to improved image quality.<br>The aluminum compensation filter provides object-specific absorption compensation for bone, tissue and air transitions. They can be attached to the multileaf collimator.                                      |

|                                                                        |                                |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                  |
|------------------------------------------------------------------------|--------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|                                                                        |                                | Weight: 0.6 kg (1.32 lbs)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                        |
| Компенсационный фильтр, плечо (не более 3 шт.).                        | Homogenizing filter shoulder   | <p>For compensation of density variations for exposures of the shoulder region. Compensation filters prevent overexposure and contribute significantly to improved image quality.</p> <p>The aluminum compensation filter provides object-specific absorption compensation for bone, tissue and air transitions. They are made of aluminum and can be attached to the multileaf collimator.</p> <p>Weight: 0.6 kg (1.32 lbs)</p>                                                                                                                                                                 |
| Держатель фильтров - 1шт. (не более 3 шт.).                            | Support for filter             | <p>For storage of compensation filters. Eight filters can be neatly stored and organized in a special holder. The holder is designed for wall installation.</p> <p>Weight: 1.3 kg (2.9 lbs)</p>                                                                                                                                                                                                                                                                                                                                                                                                  |
| Компенсационный фильтр, Т-позвоночник, L-позвоночник (не более 3 шт.). | Homogen.Filter Thor.+Lum.Sp.   | <p>To compensate for density variations in AP and lateral thoracic spine exposures, as well as lateral lumbar spine exposures. Compensating filters prevent overexposures and contribute significantly to image quality. Thoracic and lumbar spine compensating filters are used instead of graduated screens.</p> <p>The aluminum compensation filter provides object-specific absorption compensation for bone, tissue and air transitions. They can be attached to the multileaf collimator.</p> <p>Weight: 0.6 kg (1.32 lbs)</p>                                                             |
| Трехпольные шаблоны, (не более 4 шт.).                                 | Three-field templates (4 pcs.) | <p>For optical display of the IONTOMAT measurement fields on the patient by using the light localizer, so that the patient can always be properly positioned for optimal exposure. The three-field templates are inserted into the accessory rails of the primary collimator.</p> <p>Four three-field localizers for SIDs from:</p> <ul style="list-style-type: none"> <li>- 90 cm to 110 cm (35.4" to 43.3")</li> <li>- 110 cm to 130 cm (43.3" to 51.2")</li> <li>- 130 cm to 175 cm (51.2" to 68.9")</li> <li>- 175 cm to 220 cm (68.9" to 86.6")</li> </ul> <p>Weight: 0.3 kg (0.66 lbs)</p> |
| Сумка для аксессуаров.                                                 | add&swap Side Bag              | <p>This sidebag with a shoulderbag features tailored compartments for the tablet, up to two detectors and the power supply cables of your add&amp;swap upgrade kit, for easy transportation between usage points.</p>                                                                                                                                                                                                                                                                                                                                                                            |

|                                      |                                        |                                                                                                                                                                                                                             |
|--------------------------------------|----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| т – комплект<br>контроля<br>чества.  | NORMI 13 Set with holder               | Test object for quality control of X-ray installations<br>acc. DIN 6868-13. Includes a structure plate, bucky<br>mounting device, an Al absorber (directly fixed to<br>the collimator), adaption rails and a carrying case. |
| Настенное<br>зарядное<br>устройство. | Wall charger for wireless<br>detectors | Wall charger for wireless detectors                                                                                                                                                                                         |

## 6. Technical parameters;

All technical data represent typical values unless specific tolerances are given

### Patient table

| Table and tabletop           |                                                                                                                                                                                                                                                                                                               |
|------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Table height                 | 51.5 cm to 90.0 cm; total lift 38.5 cm (tabletop)                                                                                                                                                                                                                                                             |
| Tabletop dimensions          | <ul style="list-style-type: none"> <li>• Standard tabletop: 233 cm x 80 cm</li> <li>• Short tabletop (optional): 213 cm x 80 cm</li> </ul>                                                                                                                                                                    |
| Tabletop Al equivalent       | $\leq 0,7$ mm Al (with 100 kV / 3,6 mm Al HVL; IEC 60601-2-54)                                                                                                                                                                                                                                                |
| Patient weight               | Max. load 300 kg                                                                                                                                                                                                                                                                                              |
| Longitudinal travel          | <ul style="list-style-type: none"> <li>• Standard tabletop: <math>\pm 44</math> cm</li> <li>• Short tabletop (optional): <math>\pm 34</math> cm</li> </ul>                                                                                                                                                    |
| Transverse travel            | $\pm 14$ cm ( $\pm 5$ mm tolerance)                                                                                                                                                                                                                                                                           |
| Detector cover range         | $\geq 100$ cm                                                                                                                                                                                                                                                                                                 |
| Max. patient coverage        | About 190 cm with standard tabletop                                                                                                                                                                                                                                                                           |
| Grid                         | Stationary, Pb 13/92, $f_0 = 115$ cm<br>Stationary, Pb 13/40; $f_0 = 115$ cm                                                                                                                                                                                                                                  |
| Tabletop - detector distance | $\leq 73$ mm                                                                                                                                                                                                                                                                                                  |
| Maximum radiation field      | <ul style="list-style-type: none"> <li>• Maximum symmetrical radiation field at a distance of 1150 mm from the focus,</li> <li>• With table tray for MAX wi-D: 426 mm x 355 mm</li> <li>• With table tray for Core XL: 426 mm x 426 mm</li> <li>• With table tray for Core static: 426 mm x 426 mm</li> </ul> |

| Bucky wall stand                     |                                                                                                                                                                                                                                                                                                                                                                                                                             |
|--------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Travel range (central beam to floor) | From 31.5 cm to 175 cm                                                                                                                                                                                                                                                                                                                                                                                                      |
| Detector unit                        | Tiltable from $-20^\circ$ to $+90^\circ$ with detent at $0^\circ$ and $+90^\circ$                                                                                                                                                                                                                                                                                                                                           |
| Detector cover - detector distance   | <ul style="list-style-type: none"> <li>• With MAX wi-D: <math>\leq 4.2</math> cm</li> <li>• With Core XL: <math>\leq 4.2</math> cm</li> <li>• With Core static: <math>\leq 4.5</math> cm</li> <li>• With a cassette: <math>\leq 4.6</math> cm</li> </ul>                                                                                                                                                                    |
| X-ray absorption                     | $\leq 0.6$ mm Al (at 100 kV/3.6 mm Al HVL; IEC 60601-2-54)                                                                                                                                                                                                                                                                                                                                                                  |
| Anti scatter grid                    | <ul style="list-style-type: none"> <li>• Stationary grid, Pb 13/92, <math>f_0 = 140</math> cm</li> <li>• Stationary grid, Pb 13/40, <math>f_0 = 140</math> cm</li> </ul> Optional stationary grids: <ul style="list-style-type: none"> <li>• Pb 13/92, <math>f_0 = 115</math> cm; <math>f_0 = 180</math> cm; <math>f_0 = 300</math> cm</li> <li>• Pb 13/40, <math>f_0 = 115</math> cm; <math>f_0 = 180</math> cm</li> </ul> |

## X-ray tube support

| Travel range                              |                                                                                                              |        |
|-------------------------------------------|--------------------------------------------------------------------------------------------------------------|--------|
| Horizontal travel range                   | Longitudinal                                                                                                 | 352 cm |
|                                           | Transverse with 3 m trolley:                                                                                 | 215 cm |
|                                           | Transverse with 4 m trolley:                                                                                 | 350 cm |
| Vertical travel range                     | 180 cm $\pm$ 2 cm, manual or motorized                                                                       |        |
| Minimum focus-ceiling distance            | 83 cm                                                                                                        |        |
| Rotation range around the vertical axis   | -154° to + 182°, manual<br>Detents at -90°; 0°; +90°; +180°                                                  |        |
| Rotation range around the horizontal axis | $\pm$ 137° $\pm$ 2°, motorized tilting;<br>$\pm$ 120° $\pm$ 2°, manual tilting;<br>Detents at -90°; 0°; +90° |        |

| X-ray tube                                                       |                                                     |
|------------------------------------------------------------------|-----------------------------------------------------|
| Model                                                            | RAY-14S_3F                                          |
| Maximum exposure voltage (IEC 60613)                             | 150 kV                                              |
| Focal spot nominal value (IEC 60336)                             | 0.6<br>1.2                                          |
| Radiographic anode input power (IEC 60613)                       | 34 kW/80 kW                                         |
| Nominal output of anode (IEC 60613)                              | 22 kW/54 kW (50/60 Hz)<br>33 kW/78 kW (150/180 Hz)  |
| Nominal anode input power (thermal anode reference output = 0 W) | 29 kW/73 kW (50/60Hz)<br>45 kW/107 kW (150/180Hz)   |
| Anode angle                                                      | 12°                                                 |
| Maximum heat storage capacity of tube housing                    | 1,000 kJ (1,350 kHU)                                |
| Total filtration (IEC 60601-1-3)                                 | $\geq$ 2.5 mm Al/75 kV                              |
| Leakage radiation (IEC 60601-1-3)                                | (at 150kV at 1 m distance) $\leq$ 0.8 mGy/h (450 W) |

| Collimator                 |                                                                                                                                                                                                                               |
|----------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Inherent filtration        | <ul style="list-style-type: none"> <li>• Collimator AL04 II-D: 1mm Al at 70 kV</li> <li>• Collimator RFU: 1mm Al at 75 kV</li> </ul>                                                                                          |
| Full field light localizer | Very efficient high power LED technology, high energy efficiency enabling low-noise design without external cooling system; long lifetime approx. 100,000 h; timer functionality; laser line light localizer (can be covered) |
| Copper prefilter           | Without filter; 0.1 mm; 0.2 mm; 0.3 mm; motorized                                                                                                                                                                             |
| Rotation                   | $\pm$ 45°                                                                                                                                                                                                                     |
| Collimator control         | Manual or motorized                                                                                                                                                                                                           |

For Collimator AL04:

| Laser localizer |        |
|-----------------|--------|
| Class           | 2      |
| Wave length     | 650 nm |
| Power           | <1mW   |

For Collimator RFU:

| Laser localizer |        |
|-----------------|--------|
| Class           | 1      |
| Wave length     | 650 nm |
| Power           | <1mW   |

## Detectors

| MAX wi-D                                | pixium 3543 EZh                                                                                                                                                                                           |
|-----------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Detector technology                     | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology                                                                                                                        |
| Dimensions (active area)                | 34.8 cm x 42.4 cm                                                                                                                                                                                         |
| Dimensions with detector housing        | 44 cm x 46.1 cm x 1.9 cm                                                                                                                                                                                  |
| Active detector matrix                  | 2350 x 2866                                                                                                                                                                                               |
| Pixel size                              | 148 $\mu$ m                                                                                                                                                                                               |
| Semiconductor material                  | Amorphous silicon (a-Si)                                                                                                                                                                                  |
| Scintillator                            | Cesium iodide (CsI)                                                                                                                                                                                       |
| Acquisition depth                       | 16 bit                                                                                                                                                                                                    |
| Data transmission                       | < 2 s preview; < 4 s full image                                                                                                                                                                           |
| DQE in %; 2 $\mu$ Gy (RQA5) (IEC 62220) | 70% at 0.05 lp/mm<br>51% at 1 lp/mm<br>42% at 2 lp/mm<br>29% at 3 lp/mm<br>19% at Nyquist                                                                                                                 |
| MTF in % (RQA5) (IEC 62220)             | 63% at 1 lp/mm<br>35% at 2 lp/mm<br>19% at 3 lp/mm<br>12% at Nyquist                                                                                                                                      |
| Thickness                               | 19 mm                                                                                                                                                                                                     |
| Weight                                  | 3,3 kg                                                                                                                                                                                                    |
| Max. load capacity                      | <ul style="list-style-type: none"> <li>• 300 kg with patient recumbent on it (300 kg is only for functional availability rather than full performance)</li> <li>• 100 kg with patient standing</li> </ul> |
| Battery                                 | Lithium-ion, rechargeable, exchangeable                                                                                                                                                                   |
| Charging time                           | 3 h                                                                                                                                                                                                       |

|                        |                                                                                                            |
|------------------------|------------------------------------------------------------------------------------------------------------|
| Battery operation time | <ul style="list-style-type: none"> <li>Up to 1,050 images</li> <li>6.5 hours in usual operation</li> </ul> |
| Charging location      | Table detector tray, wall stand detector tray, battery charger (optional)                                  |
| Clip-on grid           | <ul style="list-style-type: none"> <li>Pb 5/85, f0 = 115 cm</li> <li>Pb 15/80, f0 = 115 cm</li> </ul>      |
| WLAN Standard          | IEEE 802.11 b/g/n, 2 x 2 mimo, WPA2/AES Encryption, EAP/TLS support                                        |

|                                   |                                                                                                                             |
|-----------------------------------|-----------------------------------------------------------------------------------------------------------------------------|
| <b>Core XL</b>                    | <b>Mars1717VS</b>                                                                                                           |
| Detector technology               | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology                                          |
| Dimensions (active area)          | 42.6 cm x 42.6 cm                                                                                                           |
| Dimensions with detector housing  | 46.1 cm x 46.1 cm x 1.57 cm                                                                                                 |
| Active detector matrix            | 3070 x 3070                                                                                                                 |
| Pixel size                        | 139 µm                                                                                                                      |
| Semiconductor material            | Amorphous silicon (a-Si)                                                                                                    |
| Scintillator                      | Cesium iodide (CsI)                                                                                                         |
| Acquisition depth                 | 16-bit                                                                                                                      |
| Data transmission                 | < 3 s preview; < 7 s full image                                                                                             |
| DQE in %; 2µGy (RQA5) (IEC 62220) | 80% at 0.05 lp/mm<br>65% at 1 lp/mm<br>53% at 2 lp/mm<br>34% at 3 lp/mm<br>21% at Nyquist                                   |
| MTF in % (RQA5) (IEC 62220)       | 64% at 1 lp/mm<br>34% at 2 lp/mm<br>18% at 3 lp/mm<br>13% at Nyquist                                                        |
| Thickness                         | 15.7 mm                                                                                                                     |
| Weight                            | 4.2 kg                                                                                                                      |
| Max. load capacity                | <ul style="list-style-type: none"> <li>150 kg with patient recumbent on it</li> <li>100 kg with patient standing</li> </ul> |
| Battery                           | Lithium-ion, rechargeable, exchangeable                                                                                     |
| Charging time                     | 4 h                                                                                                                         |
| Battery operation time            | <ul style="list-style-type: none"> <li>Up to 950 images</li> <li>7.5 hours in usual operation</li> </ul>                    |
| Charging location                 | Table detector tray, wall stand detector tray, battery charger (optional)                                                   |
| Clip-on grid                      | Pb 5/85, f0 = 115 cm                                                                                                        |
| WLAN Standard                     | IEEE 802.11 b/g/n, 2 x 2 mimo, WPA2/AES Encryption, EAP/TLS support                                                         |

|                                  |                                                                                    |
|----------------------------------|------------------------------------------------------------------------------------|
| <b>Core Static</b>               | <b>Venu1717X</b>                                                                   |
| Detector technology              | Cesium iodide scintillator coupled to TFT matrix with amorphous silicon technology |
| Dimensions (active area)         | 42.6 cm x 42.6 cm                                                                  |
| Dimensions with detector housing | 46.0 cm x 46.0 cm x 1.55 cm                                                        |
| Active detector matrix           | 3070 x 3070                                                                        |
| Pixel size                       | 139 µm                                                                             |
| Semiconductor material           | Amorphous silicon (a-Si)                                                           |

|                                   |                                                                                           |
|-----------------------------------|-------------------------------------------------------------------------------------------|
| Scintillator                      | Cesium iodide (CsI)                                                                       |
| Acquisition depth                 | 16-bit                                                                                    |
| Data transmission                 | < 3 s preview; < 5 s full image                                                           |
| DQE in %; 2μGy (RQA5) (IEC 62220) | 80% at 0.05 lp/mm<br>65% at 1 lp/mm<br>53% at 2 lp/mm<br>34% at 3 lp/mm<br>21% at Nyquist |
| MTF in % (RQA5) (IEC 62220)       | 64% at 1 lp/mm<br>34% at 2 lp/mm<br>18% at 3 lp/mm<br>13% at Nyquist                      |

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|                                          |                                                                                                                                                                                                    |
|------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Polydors RFX</b>                      |                                                                                                                                                                                                    |
| According to IEC 60601-2-7 and GB 9706.3 |                                                                                                                                                                                                    |
| High voltage form                        | Multipulse                                                                                                                                                                                         |
| Generator frequency                      | ≥ 100 kHz                                                                                                                                                                                          |
| Exposure voltage                         | 40 kV - 150 kV                                                                                                                                                                                     |
| Maximum current                          | 50 A                                                                                                                                                                                               |
| Maximum kW value in AEC mode             | 56 kW ± 5 %                                                                                                                                                                                        |
| mAs range                                | <ul style="list-style-type: none"> <li>• 0.5 mAs to 800 mAs for 55 kW/65 kW</li> <li>• 0.5 mAs to 1000 mAs for 80 kW</li> </ul>                                                                    |
| Minimum exposure time                    | <ul style="list-style-type: none"> <li>• 1-point technique: 1 ms at 100 kV</li> <li>• 2-point technique: 2 ms at 100 kV</li> <li>• 3-point technique: 20 ms at 100 kV</li> </ul>                   |
| Tolerance                                | <ul style="list-style-type: none"> <li>• kV tolerance ± 5 %</li> <li>• mAs tolerance not greater than ± (10 % + 0.2 mAs)</li> <li>• ms tolerance not greater than ± (10 % + 1 ms)</li> </ul>       |
| Nominal electric power                   | <ul style="list-style-type: none"> <li>• 55 kW : 100 kV, 55 mAs, 100 ms (550 mA )</li> <li>• 65 kW : 100 kV, 65 mAs, 100 ms (650 mA)</li> <li>• 80 kW : 100 kV, 80 mAs, 100 ms (800 mA)</li> </ul> |
| Max. electric power                      | <ul style="list-style-type: none"> <li>• 55 kW : 100 kV, 55 mAs, (550 mA )</li> <li>• 65 kW : 100 kV, 65 mAs, (650 mA)</li> <li>• 80 kW : 100 kV, 80 mAs, (800 mA)</li> </ul>                      |
| Max. tube current and voltage            | <ul style="list-style-type: none"> <li>• 55 kW: 68 kV at 800 mA</li> <li>• 65 kW: 81 kV at 800 mA</li> <li>• 80 kW: 80 kV at 1,000 mA</li> </ul>                                                   |
| mA value at nominal 150 kV               | <ul style="list-style-type: none"> <li>• 55 kW: 350 mA</li> <li>• 65 kW: 400 mA</li> <li>• 80 kW: 509 mA</li> </ul>                                                                                |


## Image station

|                        |                                                                                                                                 |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------|
| <b>Computer system</b> |                                                                                                                                 |
| Computer               | Intel® Core™ i5-6500TE Processor, 6 M Cache, up to 3.30 GHz, 1 x DDR4 8 GB memory, 2 x 500 G SATA HDD, 4 x USB 2.0, 2 x USB 3.0 |

|                                    |                                                                                                                                                                                                                                         |
|------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Operating system                   | Windows® 10 LTSB 64 bit                                                                                                                                                                                                                 |
| <b>Archiving and documentation</b> |                                                                                                                                                                                                                                         |
| System storage capacity            | 10,000 images                                                                                                                                                                                                                           |
| Background functions               | Preprocessing and data transfer is done in the background                                                                                                                                                                               |
| Archiving                          | USB, DVD or CD-R read/write in DICOM 3.0 format (or packed format)<br>CD-R storage: <ul style="list-style-type: none"> <li>• DICOM 3.0 format</li> <li>• Packed format</li> <li>• DICOM send</li> <li>• DICOM query/retrieve</li> </ul> |
| Export                             | In JPG, BMP, or TIFF format via USB/DVD/CD-R                                                                                                                                                                                            |
| Documentation                      | <ul style="list-style-type: none"> <li>• DICOM Print</li> <li>• Virtual Film Sheet</li> </ul>                                                                                                                                           |
| DICOM Worklist, MPPS               | Data exchange via HIS/RIS                                                                                                                                                                                                               |
| Paper printing                     | Printing of lists and images to a network paper printer                                                                                                                                                                                 |

## Network information

- (1) #3 NIC
- (2) #2 NIC
- (3) Imaging system (Windows 10 LTSB 2016)
- (4) #1 NIC
- (5) DICOM TCP Port 104
- (6) Switch on customer site
- (7) Any RIS or HIS or PACS

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### DAP chamber

|                                                              |       |
|--------------------------------------------------------------|-------|
| Dose area product is measured by DAP meter at the collimator |       |
| Tolerance of DAP values                                      | ± 35% |

### AEC chambers

|                                                  |      |
|--------------------------------------------------|------|
| AEC chambers attenuation equivalent of 0.7 mm Al |      |
| Tolerance of AEC values                          | ± 3% |

### SW information:

|                 |            |
|-----------------|------------|
| SW version      | VA20       |
| Date of release | 2021-03-22 |

### Patient positioning mattress:

|                            |             |
|----------------------------|-------------|
| The attenuation equivalent | < 0.5 mm Al |
| Weight                     | 2.7 kg      |
| <b>Dimensions</b>          |             |
| Width                      | 660 mm      |
| Length                     | 1980 mm     |
| Thickness                  | 25 mm       |

### Hand grip strip:

|                              |        |
|------------------------------|--------|
| The maximum permissible load | 20 kg. |
|------------------------------|--------|

### BABIX holder:

|                              |        |
|------------------------------|--------|
| The maximum permissible load | 10 kg. |
|------------------------------|--------|

### Multipurpose/Ortho Stand:

|                                                                                     |                                                      |
|-------------------------------------------------------------------------------------|------------------------------------------------------|
| Height                                                                              | 197 cm                                               |
| Depth                                                                               | 86 cm                                                |
| Width                                                                               | 82 cm (incl. brake pedals)                           |
| Weight                                                                              | 112 kg (incl. additional platform and two handgrips) |
| Max. patient weight                                                                 | 200 kg                                               |
| Max. patient size                                                                   | 190 cm                                               |
| Attenuation equivalent of the transparent protection plate acc. to DIN EN 60601-1-3 | < 1 mm Al - 100 kV 3.6 HVL                           |

### Compression Belt:

|                   |          |
|-------------------|----------|
| <b>Dimensions</b> |          |
| Width             | ≤ 100 mm |
| Length            | ≤ 400 mm |
| Height            | ≤ 150 mm |

| <b>Ambient conditions (operation)</b>             |                                                   |
|---------------------------------------------------|---------------------------------------------------|
| Temperature range                                 | +10°C to +35°C                                    |
| Relative humidity                                 | 20% to 75%, rel air humidity without condensation |
| Barometric pressure                               | 700 hPa to 1060 hPa                               |
| <b>Ambient conditions (transport and storage)</b> |                                                   |
| Temperature range                                 | -20°C to +70°C                                    |
| Relative humidity                                 | 10% to 95%, rel air humidity without condensation |
| Barometric pressure                               | 500 hPa to 1060 hPa                               |
| Aluminum equivalent                               | ≤ 0.8 mm Al                                       |

#### Cassette holder, lateral

|                          |                                                                                          |
|--------------------------|------------------------------------------------------------------------------------------|
| Maximum cassette format  | 35 cm in height, width is unrestricted, for example 35 cm x 43 cm horizontal is possible |
| Minimum cassette formats | 13 x 18 cm vertical or 18 cm x 24 cm horizontal                                          |
| Maximum load capacity    | 3 kg                                                                                     |

#### Footswitch:

| <b>Dimensions</b> |        |
|-------------------|--------|
| Width             | 157 mm |
| Length            | 310 mm |
| Height            | 35 mm  |
| <b>Weight</b>     | 2.2 kg |

#### Foot Kick Switch Front & Rear

| <b>Dimensions</b> |         |
|-------------------|---------|
| Width             | 90.5 mm |
| Length            | 715 mm  |
| Height            | 73 mm   |
| <b>Weight</b>     | 3.5 kg  |

#### Barcode reader:

| <b>Dimensions</b>     |                                       |
|-----------------------|---------------------------------------|
| Length                | 99 mm                                 |
| Width                 | 64 mm                                 |
| Height                | 165 mm                                |
| weight                | 155g                                  |
| input voltage         | 4VDC~5.5VDC                           |
| operating power       | 2.5W (500mA@5VDC)                     |
| standby power         | 0.45W (90mA@5VDC)                     |
| Interface             | USB1.1                                |
| operating environment | 0°C~50°C                              |
| humidity              | 0%~95%                                |
| Scanning mode         | Two dimensional image (1280x800pixel) |

**Three-field template, set**

- 90 cm - 110 cm
- 110 cm - 130 cm
- 130 cm - 175 cm
- 175 cm - 220 cm

**Wireless Remote Control Console**

|                                    |                                                                                                                                                                                                                                                   |
|------------------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Active range                       | Approximately 7 m around the receiver mounted to the stand.                                                                                                                                                                                       |
| Charging time                      | Approx. 8 h                                                                                                                                                                                                                                       |
| Operating time with full batteries | Approx. 20 h                                                                                                                                                                                                                                      |
| Frequency range                    | 2402 MHz - 2480 MHz                                                                                                                                                                                                                               |
| Modulation                         | GFSK                                                                                                                                                                                                                                              |
| RF Power                           | 4 dBm                                                                                                                                                                                                                                             |
| RF emissions                       | CISPR 11 Compliance: Group 2 FCC 47 CFR 15, RED, SRRC<br><br>The deployment of Wireless Remote Control in combination with MULTIX Impact C wi-D changes the classification of MULTIX Impact C wi-D from compliance group 1 to compliance group 2. |

**Tablet specification**

| Feature                           | Description                                                                                                                                                                                                                  |
|-----------------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Display type                      | Retina display                                                                                                                                                                                                               |
| Display size                      | 7.9-inch (diagonal) LED-backlit Multi-Touch display                                                                                                                                                                          |
| Screen resolution                 | 2048-by-1536 resolution at 326 pixels per inch (ppi)                                                                                                                                                                         |
| Capacity                          | 64 GB                                                                                                                                                                                                                        |
| Size (length x width x thickness) | 8.0 inches (203.2 mm) x 5.3 inches (134.8 mm) x 0.24 inches (6.1 mm)                                                                                                                                                         |
| Weight                            | 0.66 pound (300.5 grams)                                                                                                                                                                                                     |
| Chip                              | A12 Bionic chip with 64-bit architecture<br>Neural Engine<br>Embedded M12 coprocessor                                                                                                                                        |
| Camera                            | 8-megapixel camera                                                                                                                                                                                                           |
| Environmental Requirements        | Operating ambient temperature: 32° to 95° F (0° to 35° C)<br>Nonoperating temperature: -4° to 113° F (-20° to 45° C)<br>Relative humidity: 5% to 95% non-condensing<br>Operating altitude: tested up to 10,000 feet (3000 m) |

## RF transmitters for tablet


| Name            | Frequency range                                  | Bandwidth                 | Modulation type                                                                    | Effective radiated power                           |
|-----------------|--------------------------------------------------|---------------------------|------------------------------------------------------------------------------------|----------------------------------------------------|
| iPad mini A2133 | 2400-2483.5MHz,<br>5150-5350MHz,<br>5725-5850MHz | 80MHz,20MHz,8MHz,<br>3MHz | BPSK/QPSK/16QAM/<br>256QAM/DBPSK/<br>DQPSK/CCK<br><br>GFSK $\pi/4$ -DQPSK<br>8DPSK | 33dBm(Max)<br>23dBm(EIRP)(Max)<br>20dBm(EIRP)(Max) |

Siemens Healthineers hereby certifies that the Remote Interface is compliant with the following generally recognized standards covering electrical and mechanical safety:


| Standard organization | Reference number and date | Title of standard                                                                                                                                                                                                                        |
|-----------------------|---------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| IEC                   | IEC 60950-1               | Information technology equipment- Safety – Part 1: General requirements                                                                                                                                                                  |
| CISPR                 | CISPR 32                  | Electromagnetic compatibility of multimedia equipment - Emission requirements                                                                                                                                                            |
| CISPR                 | CISPR 35                  | Electromagnetic compatibility of multimedia equipment - Immunity requirements                                                                                                                                                            |
| IEC                   | IEC 62133-2               | Secondary cells and batteries containing alkaline or other non-acid electrolytes - Safety requirements for portable sealed secondary cells, and for batteries made from them, for use in portable applications - Part 2: Lithium systems |

## 7. Additional Labeling information

### MAX wi-D labeling

or

Bucky wall

**Vertic**

**Vertical module motorized**

|                     |           |
|---------------------|-----------|
| (01) 04056869272672 | GTIN      |
| (240) 11506510      | MODEL     |
| (422) 156           | ORIGIN CN |
| (21) 11005          | SERIAL    |

IVK

Siemens Shanghai Medical Equipment Ltd. 278 Zhouzhu Road  
201318 Shanghai, P.R.China

Made in P.R. China

**Tilting module Manual**

|                     |           |
|---------------------|-----------|
| (01) 04056869274263 | GTIN      |
| (240) 11506310      | MODEL     |
| (422) 156           | ORIGIN CN |
| (21) 11003          | SERIAL    |

IVK

Siemens Shanghai Medical Equipment Ltd. 278 Zhouzhu Road  
201318 Shanghai, P.R.China

Made in P.R. China

**BWS space Box 290mm**

|                     |           |
|---------------------|-----------|
| (01) 04056869274287 | GTIN      |
| (240) 11506302      | MODEL     |
| (422) 156           | ORIGIN CN |
| (21) 11005          | SERIAL    |

Siemens Shanghai Medical Equipment Ltd. 278 Zhouzhu Road  
201318 Shanghai, P.R.China

Made in P.R. China

**BWS space Box 90mm**

|                     |           |
|---------------------|-----------|
| (01) 04056869274270 | GTIN      |
| (240) 11506851      | MODEL     |
| (422) 156           | ORIGIN CN |
| (21) 11005          | SERIAL    |

Siemens Shanghai Medical Equipment Ltd. 278 Zhouzhu Road  
201318 Shanghai, P.R.China

Made in P.R. China

**Front module fixed FD**

|                     |           |
|---------------------|-----------|
| (01) 04056869272696 | GTIN      |
| (240) 11506512      | MODEL     |
| (422) 156           | ORIGIN CN |
| (21) 11005          | SERIAL    |

IVK


Siemens Shanghai Medical Equipment Ltd. 278 Zhouzhu Road  
201318 Shanghai, P.R.China

Made in P.R. China

Dose chamber labeling (CAREMAX)

Footswitch label

Symbols that can be used on the marking

|                                                                                                                                                                                       |                    |
|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|--------------------|
| GTIN                                                                                                                                                                                  | Global number GTIN |
| Model                                                                                                                                                                                 | Model              |
| ORIGIN                                                                                                                                                                                | Country of Origin  |
| <br> или SERIAL | Serial number      |
|                                                                                                    | Manufacturer       |

|                                                                                     |                                                                                                                                                                                                                                                                                                             |
|-------------------------------------------------------------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
|    | Manufacturing date                                                                                                                                                                                                                                                                                          |
|    | Do not throw it away. Hand over to a recycling center .                                                                                                                                                                                                                                                     |
|    | Compliance EC                                                                                                                                                                                                                                                                                               |
|    | Unique device identifier                                                                                                                                                                                                                                                                                    |
|    | Batch code                                                                                                                                                                                                                                                                                                  |
|    | Catalog number                                                                                                                                                                                                                                                                                              |
|    | In general, it indicates elevated, potentially dangerous levels of non-ionizing radiation. It is also used to refer to medical electrical equipment or systems that include radio frequency transmitters or purposefully use radio frequency electromagnetic energy for diagnostic or therapeutic purposes. |
|   | Attention                                                                                                                                                                                                                                                                                                   |
|  | See the operating instructions                                                                                                                                                                                                                                                                              |
|  | Temperature range                                                                                                                                                                                                                                                                                           |
|  | Use up to                                                                                                                                                                                                                                                                                                   |

### Label in Russian language

Artwork of the medical device labeling in Russian language (additionally for all the possible configurations of the MD, which are delivered on the territory of the Russian Federation, the following label in Russian is applicable):

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**8. Environment Protection Requirements;**

The system and its components comply with the following product environmental protection regulations.

| Statutory regulation | Title                                                                                                                                                                                  |
|----------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 2011/65/EU:2011      | Directive 2011/65/EU of the European Parliament and of the Council of 8 June 2011 on the restriction of the use of certain hazardous substances in electrical and electronic equipment |

07/2006/EC:2006

Regulation (EC) No 1907/2006 of the European Parliament and of the Council of 18 December 2006 concerning the Registration, Evaluation, Authorization and Restriction of Chemicals (REACH), establishing a European Chemicals Agency, amending Directive 1999/45/EC and repealing Council Regulation (EEC) No 793/93 and Commission Regulation (EC) No 1488/94 as well as Council Directive 76/769/EEC and Commission Directives 91/155/EEC, 93/67/EEC, 93/105/EC and 2000/21/EC

#### 9. Sterility (if applicable).

The system doesn't require sterilization.

#### 10. The list of International Standards / Directives / Documents applicable to the Medical Device (ISO, CE, DoC, IEC).

| Standard       | Issue | Title                                                                                                                                                                                |
|----------------|-------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| DICOM Standard | 2016  | Digital Imaging and Communications in Medicine                                                                                                                                       |
| ISO 15223-1    | 2016  | Medical devices-Symbols to be used with medical device labels, labelling and information to be supplied-Part 1:General requirements                                                  |
| EN 1041        | 2013  | Information supplied by the manufacturer with medical devices                                                                                                                        |
| EN 50581       | 2012  | Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances                                             |
| IEC 63000      | 2016  | Technical documentation for the assessment of electrical and electronic products with respect to the restriction of hazardous substances                                             |
| IEC 60601-1    | 2012  | Medical electrical equipment – General requirements for basic safety and essential performance                                                                                       |
| IEC 60601-1-2  | 2014  | Medical electrical equipment - Part 1-2: General requirements for basic safety and essential performance; Collateral standard: Electromagnetic compatibility; Requirements and tests |
| IEC 60601-1-6  | 2013  | Medical electrical equipment – Part 1-6: General requirements for basic safety and essential performance – Collateral standard: Usability                                            |
| IEC 60825-1    | 2014  | Safety of laser products - Part 1: Equipment classification and requirements                                                                                                         |
| IEC 60825-2    | 2010  | Safety of laser products - Part 2: Safety of optical fiber communication systems                                                                                                     |

|                |      |                                                                                                                                                                                                                                          |
|----------------|------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| IEC62368-1     | 2014 | Audio/video, information and communication technology equipment - Part 1: Safety requirements                                                                                                                                            |
| IEC 62040-1    | 2013 | Uninterruptible power systems (UPS) - Part 1: General and safety requirements for UPS / Combines IEC 62040-1 (2008-06) and AMD 1 (2013-01)                                                                                               |
| IEC 62304      | 2015 | Medical device software - Software life cycle processes / Combines IEC 62304 (2006-05) and AMD 1 (2015-06)                                                                                                                               |
| IEC 62353      | 2014 | Medical electrical equipment - Recurrent test and test after repair of medical electrical equipment                                                                                                                                      |
| IEC 62366      | 2014 | Medical devices - Application of usability engineering to medical devices                                                                                                                                                                |
| IEC 62366-1    | 2015 | Medical devices - Part 1: Application of usability engineering to medical devices.                                                                                                                                                       |
| IEC 62471      | 2006 | Photobiological safety of lamps and lamp systems                                                                                                                                                                                         |
| ISO 10993-1    | 2009 | Biological evaluation of medical devices - Part 1: Evaluation and testing within a risk management process                                                                                                                               |
| ISO 14155      | 2011 | Clinical investigation of medical devices for human subjects – Good clinical practice                                                                                                                                                    |
| ISO 14971      | 2007 | Medical devices - Application of risk management to medical devices                                                                                                                                                                      |
| IEC 60038      | 2009 | IEC standard voltages (Edition 7.0)                                                                                                                                                                                                      |
| IEC 80001-1    | 2010 | Application of risk management for it-networks incorporating medical devices - part 1: roles,responsibilities and activities                                                                                                             |
| ISO 17664      | 2017 | Processing of health care products - Information to be provided by the medical device manufacturer for the processing of medical devices                                                                                                 |
| IEC 60601-1-9  | 2013 | Medical electrical equipment - Part 1-9: General requirements for basic safety and essential performance - Collateral Standard: Requirements for environmentally conscious design / Combines IEC 60601-1-9 (2007-07) and AMD 1 (2013-06) |
| IEC 60601-1-3  | 2016 | Medical electrical equipment - Part 1-3: General requirements for basic safety and essential performance - Collateral standard: Radiation protection in diagnostic X-ray equipment                                                       |
| IEC 60601-2-28 | 2017 | Medical electrical equipment - Part 2: Particular requirements for the safety of X-ray source assemblies and X-ray tube assemblies for medical diagnosis                                                                                 |
| IEC 60601-2-54 | 2018 | Medical electrical equipment - Part 2-54: Particular requirements for the basic safety and essential performance of X-ray equipment for radiography and radioscopy                                                                       |

|              |      |                                                                                                                                                                                        |
|--------------|------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| EC 62220-1-1 | 2015 | Medical electrical equipment - Characteristics of digital X-ray imaging devices - Part 1-1: Determination of the detective quantum efficiency - Detectors used in radiographic imaging |
| EC 62494-1   | 2008 | Medical electrical equipment - Exposure index of digital X-ray imaging systems - Part 1: Definitions and requirements for general radiography                                          |

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#### 12. Service contacts in Russian Federation

Service contact Siemens Healthcare LLC (Customer's support center in Russian Federation):

Tel: 8 800 550 50 22

Email: [service.ru@siemens-healthineers.com](mailto:service.ru@siemens-healthineers.com)

#### 13. Warranty;

The warranty period shall begin with the day of readiness for dispatch (outgoing goods booking) for the particular product in the BA and shall be 15 months, including all logistical times. This also applies to reconditioned used systems brought back onto the market with the mark "Proven Excellence". For products listed by the BA as "replacement parts", the warranty period shall be 8 months, including all logistical times.

#### 14. Complaints;

The complaint may be submitted only for those issues which have not been the subject of the device acceptance procedure performed in accordance with the terms of the contract.

For all questions related to the device repair and maintenance please contact the manufacturer's authorized representative: Siemens Healthcare Limited Liability Company, legal and post address: 115093, Moscow, Dubininskaya str., 96 phone: 737-10-00; fax: 737-13-20, [www.siemens-healthineers.com/ru/](http://www.siemens-healthineers.com/ru/)

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