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CERTIFICATE

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China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the signature of CHEN JIANG of RAYTO LIFE AND ANALYTICAL SCIENCES CO., LTD and the seal of the said company on the annexed DOCUMENT are genuine.

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Authorized
Signature:

Zhou Siyuan

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(Date: May. 26, 2023)

证明真伪网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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
How to use this manual

Thank you for becoming the owner of the automatic coagulometer "Technology Solution 60"! Before operating the device, please read the instructions carefully.

In order to achieve optimal results, it is necessary to become well acquainted with the principle of operation of the equipment before starting clinical trials.

This manual will provide information about the installation, daily use, maintenance and other aspects of working with the automatic coagulometer "Technology Solution 60".

The differences in functions depend on the configuration version of the model. Keep all packaging materials – they are needed for storage and transportation.

For all technical questions, please contact the Authorized Representative in the Russian Federation of Laboratory Solutions LLC.

Warning: If you ignore these precautions when working with the analyzer, this may lead to distortion of the analysis results, damage to the analyzer parts and even a threat to the operator's health.

Precautions:

Note: The device is used for clinical research. A clinical diagnosis based on the results of the analysis should be carried out by a doctor in accordance with the clinical symptoms and a number of other results.

Notification

The limited liability company "Laboratory solutions" reserves the last right to make changes in the manual.

The illustrations in this manual are provided as examples and may not correspond to the actual product. The illustrations are not intended for any other purpose. Without the written permission of LLC laboratory solutions reproduction, change of this manual is forbidden.

The manufacturer is responsible for the safety, reliability and stability of the product only if the following requirements are met:

- Assembly, re-adjustment, revision and repair must be carried out by persons with special permission of LLC "Laboratory solutions»;
- the product must be operated in accordance with this manual;
- all electrical devices must comply with international standards.

Note

Specialists in the field of medical laboratory research, specially trained doctors, nurses and laboratory technicians are allowed to work with the device.

Attention

If the user is unable to ensure proper operation/maintenance, the device may fail and even pose a threat to the user's health.

The instrument must be operated under the conditions specified in this manual. Otherwise, possible undesired operation of the device, distortion of the analysis results, a damage of parts of the analyzer and injuries.

Symbols on the packaging:

Temperature limitation

Transport packaging must be kept within this temperature range.

Fragile freight

In transport packing there is fragile freight therefore it is necessary to address with it with extra care

Top

The transport package must be on the specified side up

To protect from a rain

The packed product is sensitive to moisture influence

To protect from direct sunshine

Not to roll

Symbols on the marking:

	Serial number
	Date of manufacture
	Manufacturer

	Medical device for <i>in vitro</i> diagnostics
	Compliance of the product with the basic requirements of the EU directives
	The packaging material is made of recyclable materials or can be sent for secondary processing.
	Temperature range
	Danger
	Biological hazard
	Do not touch
	Do not put your hands inside the appliance during operation.


Labeling layouts for the “Technology Solution 60” device

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


Labeling mockups for packaging

Safety precautions and recommendations

Only for in vitro diagnostics. Carefully read the following warning before use and strictly adhere to it

Attention! Before using the device, it is necessary to carefully read the safety rules and strictly observe them.

- If the appliance emits a burning smell or smoke, immediately disconnect the power supply and unplug the power plug from the outlet. If the device is used in this condition for a long time, there is a risk of fire, electric shock or injury. Contact an Authorized Representative in the Russian Federation.

- Make sure not to spill blood or reagents on the device and not to drop brackets or paper clips into the device. This may cause a short circuit. If this happens, immediately disconnect the power supply and unplug the power plug from the outlet. Then contact an Authorized representative in the Russian Federation.

- Be careful when handling samples. Be sure to use latex examination gloves; otherwise, infection is possible. If the sample gets into the eyes or a cut on the skin, rinse the contact area with plenty of water and consult a doctor immediately.

- Do not touch the internal parts when the appliance is in operation, especially moving parts, failure to do so may result in injury.

- Any amount of sample, quality control material, calibration material, etc. on the sampling probe has a potential biological hazard. Therefore, do not touch the sampling probe.

Use and disposal of reagents

- In case of contact of the reagent containing plasma with skin or clothing, rinse with plenty of clean water to prevent infection.

- If the reagent gets into your eyes, rinse them thoroughly with water and consult a doctor immediately.

- In case of accidental ingestion of the reagent, induce vomiting by drinking copious amounts of warm salt water, and immediately consult a doctor.

- In case of accidental contact with the reagent, immediately rinse the skin thoroughly with water.

- When working with samples or reagents, be sure to use gloves. After completing the work, wash your hands with a disinfectant cleaner to prevent possible infection with pathogenic and other microorganisms.

- When replacing the reagent, take precautions, wear gloves.

Voltage, connection and grounding of the power supply

- Make sure that the power supply and the grounding of the device are reliable and stable.

- Do not plug the power cord into an outlet with an AC value that does not correspond to 100-240 V. Otherwise, a fire or electric shock may occur.

- When installing, be sure to use the three-core cable supplied with the device, ensure good grounding and place the device in a place where you can freely turn off the power. Otherwise, a fire or electric shock may occur.

- Do not damage the insulation coating of the cable. Do not pull the cable or hang heavy objects with the cable. Otherwise, a short circuit may occur, which may lead to electric shock and fire.

- Turn off the power before connecting the peripheral equipment. Otherwise, electric shock or damage to the system may occur.

TABLE OF CONTENTS

Symbols on the package	3
Symbols on the marking	3
Labeling layouts for the "Technology Solution 60" device	4
Labeling mockups for packaging	7
Safety precautions and recommendations	7
Chapter 1. Introduction	12
1.1 Introduction	13
1.1.1 Name of the device	12
1.1.2 Features	13
1.2 Layout and structure of the device	13
1.2.1 Dosing system	13
1.2.2 Thermostatic system	13
1.2.3 Measuring system	14
1.2.4 Flushing system	15
1.2.5 Sample loading area	15
1.2.6 Reagent loading zone	15
1.2.7 Computer system	15
1.3 Scope of application	17
1.4 The measured parameters	17
1.5 Technical data	18
1.6. Technical characteristics of the scanner of the "Technology Solution 60" device	19
1.7. Description of medical devices intended for use in combination with the medical device "Automatic coagulometer "Technology Solution 60", not included in the delivery package.	21
Chapter 2. Installation and calibration	22
2.1 Unpacking the device	22
2.1.1 Unpacking steps	22
2.1.2 Conditions of carriage	22
2.2 Installation and operating conditions	23
2.3 Power supply requirements	23
2.4 Accession of tubes	24
2.5 Installation of the keyboard and mouse	24
2.6 Installation of the external printer	25
2.7 Calibration of the device	25
Chapter 3. Launch	26
3.1 Preparations before launch	26
3.2 Entering the system	26

3.3 Main menu	27
3.4 Functional module	28
Chapter 4. Settings of parameters	30
4.1 Settings of positions of parameters	31
4.2 Research method parameters	35
4.2.1 Coagulation method (Percentage study method)	35
4.2.2 Spectrophotometric test	36
4.2.3 Immunological analysis	36
Chapter 5. System parameters	37
5.1 Reagent parameters	37
5.2 Setting up work sheet	38
5.2.1 Choice of group of the test	39
5.2.2 Group parameters	39
5.3 The analysis with the autofilter	40
5.4 Information about the hospital	41
5.5 System parameters	42
5.6 Additional functions	43
5.6.1 Monitoring of volume of liquid of reagent	44
5.6.2 The analysis with the autofilter	44
5.7 The query log	44
5.8 Sending data	45
Chapter 6. Test	47
6.1 The position of the reagent	48
6.1.1 Setting the position of the reagent	48
6.1.2 List selection	49
6.2. Text block	49
6.2.1 Reaction unit information	49
6.2.2 Reaction unit settings	50
6.3. Sample test	51
6.3.1 Preparation for the study	51
6.3.2 Barcode scanning (if a scanner is available)	51
6.3.3 Sample collection	51
6.3.4 Explanation of the research process	52
6.3.5 Test Flow Settings	52
	53
Chapter 7. Quality control	57
7.1 Parameters Of Quality Control Materials	57
7.2 Charts	58
7.3 Quality Control Charts (Levi Jennings Charts)	58
7.4 Westgard Rules	59
7.5 Setting quality control parameters	60
7.6 Quality Control Limits	60
Chapter 8. Calibration	62

8.1 Setting of calibration parameters	63
8.2 Manual entry	65
8.3 Calibration test	66
8.3.1 Preparation for calibration	67
8.3.2 Automatic calibration dilution	67
8.3.3 Manual calibration dilution	68
8.4 Printing of calibration data	68
Chapter 9. Analytical characteristics	86
9.1 Coefficient of variation (reproducibility)	86
9.2 Sensitivity (detection limit)	86
9.3 Linearity (measuring range)	86
9.4 Influence of interferences.	86
9.5 Limitation of the method	86
9.6 Available reference materials.	86
Chapter 10. History of data	87
10.1 Inquiry according to the test	87
10.1.1 Information request according to the test	88
10.1.2 Deleting of test results	88
10.1.3 Printing of information on test parameter	88
10.2 Inquiry on the patient	89
10.2.1 Query of the patient	89
10.2.2 Change of information of the patient	90
10.2.3 Viewing of results	91
10.2.4 Deleting a patient record	91
10.2.5 Print the patient's record	91
Chapter 11. Switching off of a system	93
Chapter 12. Service of the device	94
12.1 Maintenance of the device carried out by the user	94
12.1.1 The procedure for daily maintenance	94
12.1.2 Procedure for monthly maintenance	94
12.1.3 Performing maintenance when the analyzer is idle	94
12.2 Maintenance of the device carried out by specialists	95
12.3. Replacement of the paper of the built-in printer	95
12.4 Adding a washing solution	95
12.5 Replacing the fuse	95
12.6 Software update	96
Chapter 13. System mistakes and their elimination	97
Chapter 14. Packing, storage, transportation	99
Chapter 15. Manufacturer's guarantees	100

Chapter 1 Product description

1.1 Introduction

1.1.1. Name of the device

"Automatic coagulometer "Technology Solution 60".

1.1.2. Features

1. The medical device "Automatic Coagulometer "Technology Solution 60" is a fully automatic coagulometer that collects sample data on the principle of a photoelectric sensor based on nephelometry, chromogenic method and immunoturbidimetry. It can be widely used for the diagnosis of hemorrhagic and thrombotic diseases, monitoring and therapeutic monitoring of thrombolytic and anticoagulant therapies, etc.

2. Created on the basis of the Windows XP operating system.

3. Includes the functions of automatic dilution, automatic calibration, as well as automatic display of test results, which allows you to save calibration diagrams.

4. Includes the function of determining fibrinogen, including the determination of fibrinogen by the modified Clauss method.

5. The ability to store a large amount of data: 10000 charts and search for 100,000 test results.

6. The system supports an interface for engineers to perform periodic calibration of the mechanical parts of the device.

7. Contains the auxiliary control function: built-in database of departments, doctors and system log.

8. Supports continuous addition of STAT samples during analysis.

9. HL7 medical protocol, which interacts with the laboratory's Clinical information System (LIS) in a two-way mode.

10. Supports connecting external printers that meet the following criteria.

Classification: matrix, laser, inkjet

The product includes:

1. Automatic coagulometer "Technology Solution 60" – 1 pc.

2. Instruction manual – 1 pc.

3. Device passport – 1 pc.

4. Container for distilled water – 1 pc.

5. Container for draining liquid waste – 1 pc.

6. Power cord – 1 pc.

7. Ground connection wire – 1 pc.

8. Fuse – 2 pcs.

9. Stylus – 1 pc.

10. Wire for cleaning the needle of the air intake – 1 pc.

11. Stylus holder – 1 pc.

12. Sticker for the stylus holder – 1 pc.

13. Sample cup 1.5 ml (if necessary) – no more than 100 pcs.

14. Thermal printer paper (if necessary) – no more than 2 pcs.

15. Adapter for samples (if necessary) – no more than 1 pc.

16. Adapter for reagents (if necessary) – no more than 3 pcs.
17. Reagent bottle □23 mm (if necessary) – no more than 2 pcs.
18. Reagent bottle □35 mm (if necessary) – no more than 2 pcs.
19. Reagent bottle □28 mm (if necessary) – no more than 2 pcs.
20. Magnetic stirrer (if necessary) – no more than 2 pcs.
21. Sticker for reagents (if necessary) – no more than 1 pc.
22. Sample reader (if necessary) – 1 pc.
23. Reagent reader (if necessary) – 1 pc.

NOTE: The use of a barcode scanner is optional. To resolve the installation issue, please contact an Authorized representative of the manufacturer in the Russian Federation.

1.2 Device and accessories

The device has a number of fully automatic functions, such as adding reagent and sample, counting, displaying and printing test results, etc. The instrument mainly consists of a sample monitoring system, a constant temperature maintenance system, a measuring system, a computer system and a flushing system.


Рисунок 1-1. Front panel view

1.2.1 Dosing system

The dosing system performs automatic addition of samples and reagents and performs the following functions:

1. Determination of the liquid level. If there is not enough reagent or sample, the device will issue a corresponding message and thereby notify the user.
2. Preheating function, which allows the reagent to be heated to the required temperature in 3-5 seconds.

1.2.2 Thermostatic system

The thermostatic system includes cooling of reagents, preheating of the reagent, incubation of the reaction liquid and maintaining the temperature regime during testing.

1. Incubation function: contains eight incubation cells in the detection zone with constant temperature maintenance of 37 °C. Control over the constant temperature of sample

positions: contains a total of seven sample positions in the detection zone with a constant temperature of 37 °C.

2. Cooling: the temperature of all reagent positions is $\leq 15^{\circ}\text{C}$.

3. Reagent preheating: the temperature of the reagent probe is maintained at 37°C. The reagent is preheated in 3-5 seconds.

1.2.3 Measuring system

The measurement system turns on the light flux, collects signals and preprocesses them during the reaction. The measuring system includes three main testing methods: coagulation, immunoturbidimetric and chromogenic methods.

1. Coagulation method: five test positions are used to perform it: 9, 10, 11, 12 and 13, as shown in Figure 1-2. The wavelength is 660 nm.

2. Chromogenic method: one position is involved. The test wavelength is 405 nm, marked in Figure 1-2 with the number 14.

3. Immunoturbidimetric method: one position is involved. The test wavelength is 575 nm, marked in Figure 1-2 with the number 15.

4. 8 incubation positions, marked in Figure 1-2 with the numbers 1-8.

5. The position of the discharge of spent cuvettes is marked in Figure 1-2 with the number 16.

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To obtain correct analysis results and correct operation of the device, it is necessary to use a medical device in accordance with RU No. RZN 2020/10843 of June 17, 2020 "Disposable cuvette for the "Technology Solution" line of coagulometers", in versions:

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 μl) – 1000 pcs. in the package.
- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 μl) – 1500 pcs. in the package.

Technical characteristics of the cuvette:

- The total height of the cuvette is 30.06 ± 0.15 mm;
- Diameter of the landing part of the cuvette - 10.00 ± 0.1 mm;
- External diameter of the bottom of the cuvette - 7.70 ± 0.1 mm;
- Distance of the measurement area from the mating plane - 21.40 ± 0.1 mm;
- External diameter of the measuring part of the cuvette - 7.75 ± 0.1 mm;
- Wall thickness in the measuring part of the cuvette is not more than - 1.08 ± 0.03 mm;
- The maximum deviation of the cell wall thickness in the measuring part is $0 \dots + 0.01$ mm.
- The displacement radius is 1.60 ± 0.03 mm.

Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 has a special marking on the upper face of each cuvette:

In order to prevent the use of counterfeit medical products, it is prohibited to use cuvettes K2 without the specified marking on the line of coagulometers "Technology Solution". If the cuvette is used without the appropriate marking, the authorized representative in the Russian Federation of Laboratory Solutions LLC disclaims any responsibility for the quality of the results and for the operability of the device, and also has the right to refuse replacement under warranty of the failed device or any part of it.

1.2.4 Flushing system

The flushing system includes.

1. Flushing tank: The left hole is used to flush the inside of the sampler, and the right hole is used to flush the outside of the sampler.
2. Washing tank: Used for storing distilled or ionized water.
3. Waste container: used to collect liquid waste.
4. As a detergent and disinfectant reagent for cleaning the coagulometer dosing system, a medical device must be used in accordance with registration certificate No. RZN 2019/9048 dated May 04, 2021 "Washing solution according to technical specifications 21.20.23-001-34065548-2019".

1.2.5 Sample loading area

The sample loading area is represented by three interchangeable tripods of 9 positions, and any free position of the sample can be used to install an urgent sample. Standard test tubes, blood sampling tubes (vacutainers) and a 1.5 ml sample container can be used for the study. Blood sampling tubes (vacutainers) can have an outer diameter of 10 mm ~ 16 mm and no more than 100 mm in height.

1.2.6 Reagent loading zone

The reagent loading zone is represented by a removable tripod at 23 positions, where the first three positions are reserved for diluents and washing solution. As well as a system

for maintaining the temperature below 15 ° C. The temperature stabilization process takes 30 seconds.

In the reagent tripod, containers with a volume from 1 ml to 10 ml with an outer diameter from 14 mm to 36 mm can be installed. The height of the reagent container should not exceed 80 mm, otherwise damage may occur

1.2.7 Computer system

The system coordinates various actions and controls the operation of the device, including the interaction of man and equipment, displays information and performs printing.

1.3 Scope of application

The analyzer uses the principle of optical detection, which allows performing clotting, chromogenic and immunoturbidimetric techniques for analyzing the time and activity of blood clotting factors, as well as analyzing the activity of anticoagulant and fibrinolytic systems.

The presence of detection methods (clotting, chromogenic, immunoturbidimetry) allows you to use any test system based on these methods on the device, according to the parameters specified in the instructions of the test systems used. Thereby significantly expanding the types of tests specified in clause 1.4. of this manual.

Human blood plasma stabilized with an anticoagulant (sodium citrate) is used for analysis. Any other material is considered unintended.

1.4 The measured parameters

Measured parameters are shown in Table 1-1.

Table 1-1

Name	Measured Parameter	Analysis Method	Units
Prothrombin time	PT	Coagulation	s (second)
Activated partial thromboplastin time	APTT	Coagulation	s (second)
Fibrinogen (standard Clauss method)	FIB	Coagulation	g / l
Fibrinogen (modified Clauss method)	FIB MMC	Coagulation	g / l
Thrombin time	TT	Coagulation	s (second)
Lupus Anticoagulant LA1	LA1	Coagulation	s (second)
Lupus Anticoagulant LA2	LA2	Coagulation	s (second)
Coagulation factor	II, V, VII, VIII, IX, X, XI, XII	Coagulation	%
Coagulation factor	XIII	Chromogenic method	%
Antithrombin III	AT3	Chromogenic method	%
Plasminogen	Plg	Chromogenic method	%
Protein C	PCch	Chromogenic method	%

Protein C	PCco	Coagulation	%
Protein S	PSch	Chromogenic method	%
Protein S	PSco	Coagulation	%
Heparin	Hep	Chromogenic method	IU / ml
DOAC	Hep2	Chromogenic method	ng / ml
α 2 antiplasmin	α 2 AP	Chromogenic method	%
Fibrinogen degradation products	FDP	Immunoturbidimetry	mkg / ml
D-Diimer	D-dimer	Immunoturbidimetry	ng / ml
Antigen factor Von Willebrand	vWF Ag	Immunoturbidimetry	%
Fibrin-monomer	FM	Immunoturbidimetry	mkg / ml

1.5 Technical data

Table 1-2

Testing method	Coagulation process (optical scattering method), chromogenic method and immunoturbidimetry (optical transmission method) Coagulation process (optical scattering method), chromogenic method and immunoturbidimetry (optical transmission method)
Wavelengths	Coagulation process: 660 nm Chromogenic method: 405 nm Immunoturbidimetry: 575 nm
Measuring Positions	7 positions
Incubation positions	8 positions
Discharge positions for ditches:	1 position
Performance:	60 tests / hour (PT test)
Sample positions:	Sample stand with 27 test positions
Reagent positions:	23 cooled positions for reagents
Maximum load of measuring cuvettes:	72 measuring cuvettes
Nutrition:	a.c. 100-240 V, 50 Hz / 60 Hz
Power:	400 watts
Terms of Use:	10 °C ~ 30 °C, relative humidity ≤70%
Storage conditions:	- 20 °C 55 °C, relative humidity ≤85%, Pressure: 50 kPa ~ 106 kPa, in a well ventilated

	room
Interface:	RS-232, WLAN port, USB port (2), VGA output, PS / 2 connectors for keyboard and mouse
The maximum allowable time for setting the operating mode, calculated from the moment the analyzer is turned on	up to 10 minutes
Safety class for software acc. To EN 62304	A
Electrical safety according to EN 61010-1	- Overvoltage category I. - Pollution Degree 1. - Protection class I.
Group of dependence on perceived mechanical influences according to GOST R 50444-2020	2 - wearable, portable and mobile, not intended for use when carrying and moving within a stationary room
The degree of protection of the analyzer from dangerous penetration of water or solid particles - IP	IP 00

1.6 Description of the composition of the automatic coagulometer "Technology Solution 60"

№ п/п	Description	Assignment	Technical specifications
1.	Capacity for distilled water	It is intended for a distilled water	Volume: 4 liters. Overall dimensions: Height 325 mm $\pm$ 1 mm; Width 160 mm $\pm$ 1 mm; Depth 115 mm $\pm$ 1 mm
2.	Tank for discharging liquid waste	It is intended for collecting the fulfilled liquid	Volume: 4 liters. Overall dimensions: Height 325 mm $\pm$ 1 mm; Width 160 mm $\pm$ 1 mm; Depth 115 mm $\pm$ 1 mm.
3.	1.5 ml sample cup (if necessary)	Are used for samples or scarce plasmas during test by the device	Volume: 1.5 ml. Overall dimensions: Diameter: 15mm $\pm$ 1mm; Height: 25 mm $\pm$ 1 mm.
4.	Power cord	Provides electricity feed-in to the device	Length 2000 mm $\pm$ 10 mm. Voltage: up to 500 V. Diameter: 7 mm $\pm$ 0.5 mm.
5.	Ground wire	Provides protective grounding of the device	Voltage: up to 600 V. Length: 5200 mm $\pm$ 10 mm. Diameter: 4 mm $\pm$ 0.5 mm.
6.	Thermal printer paper	It is intended for display of	The maximum diameter

	(if necessary)	results of testing on paper	of the tape: 50 mm. Tape width: 57 mm. The paper density is 80 g/m ² .
7.	Sample Adapter (if required)	It is intended for fixing ditches of samples in a support	Diameter: 24 mm ± 1 mm, Height: 29 mm ± 1 mm.
8.	Adapter for reagents (if necessary)	It is intended for fixing of a bottle of reagents in a support	Diameter: 24 mm ± 1 mm, Height: 29 mm ± 1 mm.
9.	Fuse	Performs protective function, protects an electrical circuit and its elements from overheating and ignition at course of current of high force.	Diameter x Length: 5x20 mm. Current strength: 5 A. Voltage: 250 V.
10.	Wire for cleaning the needle of the sampler	It is used for a clearing of a needle of a probozabornik in case of contamination	Length: 445 mm. Diameter: 0.5 mm.
11.	Stylus	It is used for management of the touch interface of the coagulometer	Type: capacitive. Length: 100 mm.
12.	Stylus Holder	It is necessary for convenient fixing of the stylus	Length: 65 mm. Width: 15 mm.
13.	Sticker for stylus holder	It is necessary for fixing of fastening of the stylus on the surface of the device	Double-sided adhesive tape.
14.	Reagent bottle Ø23 mm (if necessary)	It is intended for reagent storage	Volume: 25 ml. Overall dimensions: Diameter: 23mm ± 1mm; Height: 45 mm ± 1 mm.
15.	Reagent bottle Ø35 mm (if necessary)	It is intended for reagent storage	Volume: 30 ml. Overall dimensions: Diameter: 35mm ± 1mm; Height: 60 mm ± 1 mm.
16.	Reagent bottle Ø28 mm (if necessary)	It is intended for reagent storage	Volume: 20 ml. Overall dimensions: Diameter: 28mm ± 1mm; Height: 50 mm ± 1 mm.
17.	Magnetic stirrer (if necessary)	It is intended for hashing of reagents	Length: 15 mm.
18.	Reagent Label (if required)	It is used for designation of installation sites of reagents in a support	Length: 13 mm. Width: 5 mm.
19.	Sample reader (if necessary)	It is intended for reading of a barcode of samples	Specifications: • Wavelength: 650 nm. • Power: up to 1 MW. • Laser class: 2. • Maximum distance to the barcode: 90 cm. • Reading speed of 104 ±

			12 barcodes per second. • Power consumption: 95 MW. • System connection: hardware.
20.	Reagent Reader (if necessary)	It is intended for reading of a barcode of reagents	Specifications: • Wavelength: 650 nm. • Power: up to 1 MW. • Laser class: 2. • Maximum distance to the barcode: 90 cm. • Reading speed: 104 ± 12 barcodes per second. • Power consumption: 95 MW. • System connection: USB. • Length

1.7. Description of medical devices intended for use in combination with the medical device "Automatic coagulometer "Technology Solution 60", not included in the delivery package.

Name	Number and date of issue of the registration certificate	Appointment
A set of reagents for the determination of activated partial thromboplastin time for the line of coagulometers "Technology Solution" (TS-APTV) according to TU 21.20.23-082-42349142-2020	RZN 2020/12933 of December 14, 2020	The medical device "TS-APTV" is designed to perform a basic technique for the study of the hemostasis system - the determination of activated partial (partial) thromboplastin time (APTV / APTT) on the line of semi-automatic and automatic coagulometers "Technology Solution". APTT is determined to evaluate the internal pathway of coagulation activation in blood plasma.
A set of reagents for determining the prothrombin time for the line of coagulometers "Technology Solution" (TS-Thromboplastin) according to TU 21.20.23-083-42349142-2020	RZN 2020/12930 of December 15, 2020	The medical device "TS-Thromboplastin" is designed to determine the prothrombin clotting time (PV) in venous blood plasma, calculate the prothrombin ratio (PO), the prothrombin index (PTI), the international normalized ratio

		(INR) and the prothrombin index by Quick on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of prothrombin time is carried out when monitoring the treatment with anticoagulants of non-direct action, as well as to detect dysfunction of the external coagulation pathway.
A set of reagents for determining thrombin time for the line of coagulometers "Technology Solution" (TS-Thrombin) on the TU 21.20.23-084-42349142-2020	RZN 2020/12659 of November 24, 2020	The medical device "TS-Thrombin" is designed to determine the thrombin time on the line of semi-automatic and automatic coagulometers "Technology Solution". Determination of thrombin time is carried out to assess the final stage of blood plasma coagulation during prophylactic or therapeutic administration of anticoagulants, in connection with arterial or venous thrombosis (myocardial infarction, deep vein thrombosis, pulmonary embolism, etc.), as well as in the detection of congenital or acquired (cirrhosis of the liver) dysfibrinogenemia.
A set of reagents for determining the concentration of fibrinogen for the line of coagulometers "Technology Solution" (TS-Fibrinogen) on the TU 21.20.23-085-42349142-2020	RU No. RZN 2020/12640 of November 30, 2020	The medical device "TS-Fibrinogen" is intended for quantitative determination of fibrinogen in blood plasma without preliminary dilution of the studied plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". Determination of fibrinogen concentration is carried out with bleeding of various etiologies (during pregnancy and childbirth, against the background of fibrinolytic therapy, etc.), with inflammatory, autoimmune, infectious and septic diseases

		(pyelonephritis, pneumonia, peritonitis, sepsis, systemic lupus erythematosus, etc.).
A set of calibrators for determining the concentration of fibrinogen for the line of coagulometers "Technology Solution" (TS-Fibrinogen calibrator) according to TU 21.20.23-086-42349142-2020	RZN 2020/12691 of November 30, 2020	The medical device "TS-Fibrinogen-calibrator" is designed to obtain calibration values of clotting time when determining the concentration of fibrinogen in blood plasma by the modified Clauss method without preliminary dilution of the plasma under study on the line of semi-automatic and automatic coagulometers "Technology Solution". The medical device "TS-Fibrinogen-calibrator" is intended only for the set "TS-Fibrinogen".
A set of reagents for the determination of antithrombin for the line of coagulometers "Technology Solution" (TS-Antithrombin) on the TU 21.20.23-087-42349142-2020	RZN 2020/12686 of November 30, 2020	The medical device "TS-Antithrombin" is designed to determine the activity (as a percentage of the norm) of the physiological anticoagulant – antithrombin (AT) on the line of coagulometers "Technology Solution". Determination of AT is carried out to diagnose congenital and acquired deficiency of this anticoagulant (DIC syndrome, various types of shock, sepsis, obstetric pathology, severe trauma, etc.), as well as with prolonged or high-dose use of heparins.
A set of reagents for the quantitative determination of D-dimer in blood plasma for the line of coagulometers "Technology Solution" (TS-D-dimer) on the TU 21.20.23-088-42349142-2020	RZN 2020/12270 of November 14, 2020	The medical device "TS-D-dimer" is intended for the quantitative determination of D-dimer in human blood plasma by the immunoturbidimetric method at a wavelength of 570-800 nm on the line of coagulometers "Technology Solution". D-dimer is one of the most important

		markers of activation of blood clotting and fibrinolysis, since it is formed in the bloodstream as a result of the formation and subsequent lysis of stabilized fibrin by plasmin. A high level of D-dimer will be registered in patients with thrombosis, thromboembolism, DIC syndrome of various genesis.
TS-D-dimer control	RZN 2020/12571 of November 17, 2020	The medical device "TS-D-dimer control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of D-dimer in blood plasma for the line of coagulometers "Technology Solution" (TS-D-dimer)" on the line of coagulometers "Technology Solution".
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS-calibrator) on the TU 21.20.23-090-42349142-2020	RZN 2020/12354 of October 28, 2020	The medical device "TS-calibrator" is used to obtain calibration values and build calibration curves for the study of the hemostasis system on the line of semi-automatic and automatic coagulometers "Technology Solution". The use of calibration plasma frees in most cases from the need to receive fresh normal blood plasma from healthy people. Calibration plasma "TS-calibrator" is certified for constructing calibration curves and obtaining calibration values when determining the following indicators: - APTV/ACTV;

		- prothrombin time; - prothrombin index by Quick; - thrombin time; - antithrombin; - plasminogen; - protein C; - Willebrand factor antigen level (vWF:Ag); - the level of protein S (free) antigen.
A set of control plasmas for the line of coagulometers "Technology Solution" (TS-control) according to TU 21.20.23-091-42349142-2020	RZN 2020/12097 of September 30, 2020	The medical device "TS-control" is designed for in-laboratory quality control (Quality Control) in the study of the hemostasis system on the line of semi-automatic and automatic coagulometers "Technology Solution". The set is certified in the normal and pathological ranges according to the following indicators: - APTV/ACTV; - prothrombin time; - international normalized attitude (INR); - prothrombin index by Quick; - thrombin time; - - fibrinogen (modified by Klaus method); - antithrombin; - plasminogen; - protein C;

		- factor VIII; - factor IX; - Willebrand factor antigen level (vWF:Ag); - the level of protein S (free) antigen.
A set of reagents for determining the activity of factor VIII for the line of coagulometers "Technology Solution" (TS-Factor VIII) according to TU 21.20.23-092-42349142-2020	RZN 2020/13132 of January 12, 2021	The medical device "TS-Factor VIII" is designed to determine the activity of factor VIII in blood plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of factor VIII is used for the diagnosis of hemophilia A, for monitoring the replacement therapy of patients with hemophilia A with factor VIII compounds, as well as for the diagnosis of thrombophilia caused by an increased level of factor VIII.
A set of reagents for determining the activity of factor IX for the line of coagulometers "Technology Solution" (TS-Factor IX) according to TU 21.20.23-093-42349142-2020	RZN 2021/13310 of January 28, 2021	The medical device "TS-Factor IX" is designed to determine the activity of factor IX in blood plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of factor IX is used to diagnose hemophilia B, to control the replacement therapy of patients with hemophilia B with factor IX concentrates.
A set of reagents for the determination of plasminogen for the line of coagulometers "Technology Solution" (TS-Plasminogen) on the TU	RZN 2021/13401 of February 08, 2021	The medical device "TS-Plasminogen" is designed to determine the activity (as a percentage of the norm) of the main component of the fibrinolytic system – plasminogen

21.20.23-095-42349142-2020		on the line of coagulometers "Technology Solution". Determination of plasminogen is used in the diagnosis of DIC syndrome and thrombophilia, control of treatment with drugs that activate fibrinolysis.
A set of reagents for the determination of protein C for the line of coagulometers "Technology Solution" (TS-Protein C) according to TU 21.20.23-096-42349142-2020	RZN 2021/13611 of March 02, 2021	The medical device "TS-Protein C" is designed to determine the activity (as a percentage of the norm) of protein C in blood plasma on the line of coagulometers "Technology Solution". Protein C is one of the most important physiological anticoagulants synthesized in the liver with the participation of vitamin K. In its active form, it destroys coagulation factors VIIIa and Va. Hereditary protein C deficiency leads to thrombosis and thromboembolism. Newborns with protein C deficiency may develop fulminant purpura.
A set of reagents for the determination of lupus anticoagulant for the line of coagulometers "Technology Solution" (TS-Lupus) according to TU 21.20.23-094-42349142-2020	RZN 2021/14378 of May 20, 2021	The medical device "TS-Lupus" is intended for the determination of lupus anticoagulant (VA) on the line of semi-automatic and automatic coagulometers "Technology Solution" in sodium citrate-stabilized blood plasma. In blood plasma, VA binds to negatively charged phospholipids and protein-phospholipid complexes and inhibits the activation and interaction of plasma coagulation factors. These disorders are most clearly detected in phospholipid-dependent coagulation tests. The presence of VA in plasma is accompanied by recurrent

		thrombosis of veins and arteries, disorders of cerebral circulation (headaches, fainting, dynamic disorders of cerebral circulation, paresis, episynodrome, visual disturbances, etc.), fetoplacental insufficiency, habitual miscarriage (miscarriages, fetal death of the fetus), thrombocytopenia, less often – microcirculatory bleeding, polyallergia, others immune disorders, a tendency to develop DIC syndrome.
A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa) on the TU 21.20.23-102-42349142-2021	RZN 2021/15119 of August 18, 2021	The medical device "TS-Anti-Xa" is intended for quantitative determination of the anti-Xa activity of direct anticoagulants in blood plasma by chromogenic method on the line of coagulometers "Technology Solution" when monitoring the therapy with anticoagulants, such as unfractionated and low molecular weight heparins, direct oral anticoagulants (DOAC) (rivaroxaban, apixaban, fondaparinux and others).
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS-Heparin calibrator) on the TU 21.20.23-103-42349142-2021	RZN 2021/15136 of August 20, 2021	The medical device "TS-Heparin calibrator" is designed to obtain calibration values for determining the anti-Xa activity of direct anticoagulants (unfractionated and low molecular weight heparins) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution". The medical device "TS-Heparin calibrator" is intended only for the "TS-Anti-Xa" kit.
A set of control plasmas for the line of coagulometers	RZN 2021/15104 of	The medical device "TS-Heparin-control" is intended for in-

"Technology Solution" (TS-Heparin control) according to TU 21.20.23-105-42349142-2021	August 16, 2021	laboratory quality control of a medical device "A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa)" when determining the anti-Xa activity of direct anticoagulants (unfractionated heparin (NG) and low molecular weight heparins (NMH)) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution".
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS-DOAC -calibrator) on the OTHER 21.20.23-104-42349142-2021	RZN 2021/15411 of September 23, 2021	The medical device "TS- DOAC -calibrator" is designed to obtain calibration values for determining the anti-Xa activity of direct oral anticoagulants (DOAC) (apixaban, rivaroxaban, fondaparinux and others) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution". The medical device "TS- DOAC -calibrator" is intended only for the set "TS-Anti-Xa".
A set of control plasmas for the line of coagulometers "Technology Solution" (TS-DOAC -control) according to TU 21.20.23-106-42349142-2021	RZN 2021/15147 of August 20, 2021	The medical device "TS- DOAC -control" is intended for in-laboratory quality control of a medical device "A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa)" when determining the anti-Xa activity of direct oral anticoagulants (DOAC) (apixaban, rivaroxaban, fondaparinux and others) in blood plasma by chromogenic method on the line of coagulometers

		"Technology Solution".
A set of reagents for the determination of the Willebrand factor antigen for the line of automatic coagulometers "Technology Solution" (TS-Willebrand Ag) according to TU 21.20.23-107-42349142-2021	RZN 2021/15490 of October 04, 2021	The medical device "TS-Willebrand Ag" is intended for the quantitative determination of the Willebrand factor antigen (vWF:Ag) in blood plasma. Willebrand factor (vWF) is a large multimeric glycoprotein whose main functions are participation in platelet adhesion by binding it to glycoprotein Ib, as well as protection of antihemophilic globulin from proteolysis. Definition of vWF:Ag is one of the most popular methods of diagnosing hemorrhagic disease, which is based on a deficiency or functional inferiority of vWF (Willebrand's disease, Willebrand syndrome). In addition, it is advisable to use this method to monitor the effectiveness of the treatment of this pathology with Willebrand factor concentrates.
A set of reagents for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer) on the TU 21.20.23-113-42349142-2021	RZN 2022/16545 of February 15, 2022	The medical device "TS-Fibrin-monomer" is intended for the quantitative determination of fibrin monomers (FM) in human blood plasma either in free form or as part of soluble fibrin-monomer complexes (RFMC) by immunoturbidimetry on coagulometers "Technology Solution". The definition of FM is important in various variants of intravascular blood clotting (thrombosis of the main vessels, PE, DIC syndrome), including the use of anticoagulants.
A set of control plasmas for quantitative determination of FM on coagulometers	RZN 2022/19288 of December 30, 2022	The medical device "TS-Fibrin-monomer control" is intended for in-laboratory quality control of a

"Technology Solution" (TS-Fibrin-monomer control) according to TU 21.20.23-115-42349142-2021		medical device "A set of reagents for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer)". The medical device "TS-Fibrin-monomer control" is intended only for the set "TS-Fibrin-monomer".
A set of calibration plasmas for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer calibrator) according to TU 21.20.23-117-42349142-2021	RZN 2022/19253 of December 27 , 2022	The medical device "TS-Fibrin-monomer calibrator" is designed to obtain calibration values for the quantitative determination of fibrin-monomer (FM) with a set of "TS-Fibrin-monomer" on coagulometers "Technology Solution".
A set of reagents for quantitative determination of FDP on coagulometers "Technology Solution" (TS- FDP) according to TU 21.20.23-114-42349142-2021	RZN 2022/19267 of December 28 , 2022	The medical device "TS- FDP " is designed for the quantitative determination of fibrinogen and fibrin degradation products (FDP) in human blood plasma by immunoturbidimetry on coagulometers "Technology Solution". FDP s are a marker of activation of blood clotting and fibrinolysis, since they are formed in the bloodstream as a result of plasmin lysis of fibrinogen and fibrin. A high level of FDP will be registered in patients with thrombosis, thromboembolism, DIC syndrome of various genesis
A set of calibration plasmas for quantitative determination of PDF on coagulometers "Technology Solution" (TS-FDP calibrator) on the TU 21.20.23-118-42349142-2021	RZN 2022/19286 of December 30, 2022	The medical device "TS- FDP control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of FDP on coagulometers "Technology Solution" (TS-FDP)". The medical device "TS-FDP control" is intended only for

		the set "TS-FDP ".
A set of control plasmas for quantitative determination of FDP on coagulometers "Technology Solution" (TS-FDP control) according to TU 21.20.23-116-42349142-2021	RZN 2022/19286 of December 30, 2022	The medical device "TS- FDP control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of FDP on coagulometers "Technology Solution" (TS-FDP)". The medical device "TS-FDP control" is intended only for the set "TS- FDP ".
Disposable cuvette for the line of coagulometers "Technology Solution", in versions: - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 µl) – 1000 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 µl) – 1500 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K3 (1 ml) - 1000 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K4 (1.5 ml) - 500 pcs. in the package. - Disposable cuvette for	RZN 2020/10843 of June 17, 2020	The medical device "Disposable cuvette for the line of coagulometers "Technology Solution"" is intended for placing blood plasma samples and reagents, mixing them and examining blood plasma samples using coagulometric analysis methods on the line of coagulometers "Technology Solution": "Technology Solution 4", "Technology Solution 60", "Technology Solution 190", "Technology Solution 400".

coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K 5 (2 ml) - 400 pcs. in the package.		
"Washing solution according to TU 21.20.23-001-34065548-2019 "	RZN 2019/9048 of May 04 , 2021	The medical device "Washing solution" is used as a detergent and disinfectant reagent for cleaning the dosing systems of the line of coagulometers "Technology Solution" ("Technology Solution 60", "Technology Solution 190", "Technology Solution 400"), coagulometers "AK" and "Insula 500".

When new sets of the Technology Solution line appear, users will be informed by an authorized representative in the Russian Federation of Laboratory Solutions LLC through an Official Notification, or on the official website www.labsolution.ru , www.лаб-решение.RF.

Chapter 2 Installation and calibration

2.1 Unpacking the device

2.1.1 Unpacking steps

Unpack the appliance and remove the material for transportation. Keep all packaging materials and materials for transportation in case of possible transportation of the device.

The box with the device should stand straight, and the arrows should be pointing upwards.

Use the tools to open the side walls of the drawer, take out the accessories and check the package according to the packing list. In case of a shortage, contact the distributor or representatives of the Limited Liability Company "Laboratory Solutions" immediately. Check the contents, which must correspond to the list specified in clause 1.1

Remove the packaging material from the box and carefully lift the device by special recesses on both sides with the help of at least two people.

Attention: Accessories must correspond to the list of equipment in the device passport. If any component is missing, damaged or does not match the configuration in the device passport, contact an Authorized Representative in the Russian Federation. Please store the packaging properly for further packaging and transportation. The ambient transport environment for the device: temperature $-20\text{ }^{\circ}\text{C} \sim 55\text{ }^{\circ}\text{C}$, relative humidity up to 90%, barometric pressure $50\text{ kPa} \sim 106\text{ kPa}$. During transportation, pay attention to the icons "Top" "Fragile" "Protect from rain and moisture" "Limit of layers on the package" and transport the device according to the requirements.

2.1.2 Conditions of carriage

- 1) Before transporting the instrument, make sure that it is reset, the sampling needle is in the raised position.
- 2) Under normal conditions, trolleys and other vehicles can be used for transport over short distances.
- 3) During the entire transportation of the instrument, pay special attention to the safety of the display and sampling needle.
- 4) During transportation it is recommended to keep the device in an upright position, do not tilt.
- 5) Avoid vibrations if possible. Before the subsequent operation of the device, it is recommended to check and calibrate.

2.2 Installation and operating conditions

The automatic coagulometer Technology Solution 60 must be installed by qualified personnel. When it is necessary to change the installation location of the device, also contact qualified specialists.

It should be noted that if malfunctions have occurred, as a result of the movement of the system carried out by the user, this leads to the loss of all warranty obligations, even if the warranty period has not expired.

In order to ensure the correct operation of the device, it is necessary to place it on a work surface that meets the following requirements:

- Away from direct sunlight;
- In a dust-free place;
- Absence of strong electromagnetic radiation;
- Do not use this device at low relative humidity, especially if synthetic materials (synthetic clothing, carpets, etc.) are present near it, may lead to erroneous results due to the influence of electrostatic discharges.
- A grounded power outlet will be located in the immediate vicinity of the appliance;
- Leave enough space for work and maintenance. Creating conditions for heat exchange of the working device, leave approximately 50 cm from the walls to the back panel, sides, top of the device;
- The room should be well ventilated;
- Avoid excessive humidity and high temperatures;
- Avoid vibration and shaking.
- The device must be installed on a flat working surface, without inclination and irregularities.

The manufacturer is responsible for providing the consumer or customer with information about the electromagnetic compatibility of the equipment.

The consumer is responsible for maintaining an electromagnetic environment for the equipment that ensures compatibility, in which the equipment must function in accordance with its purpose.

Attention:

This device operates at an ambient temperature of 15 °C ~ 30 °C and humidity ≤ 70%.

After installation, avoid moving the appliance frequently. Use a stable trolley to move it. The angle of inclination should not exceed 15°.

Installation and relocation must be carried out exclusively by qualified specialists.

2.3 Power supply requirements


a.c. 100-240 V; 2) 50 Hz / 60 Hz; 3) 400 watts

Warning:

- 1) *The AC power supply must be well grounded (ground voltage < 5V).*
 - 2) *The power supply must be stable when combined with high-power electrical appliances.*
 - 3) *If smoke, a sharp smell or unusual noise emanating from the device is detected, immediately turn off the power and contact the service center of an Authorized representative in the Russian Federation.*
 - 4) *To disconnect the power cord, gently pull the plug, holding directly to the cord is highly discouraged.*
 - 5) *Do not open the back cover, side covers and panels of the switched-on appliance.*
- 2.4 Accession of tubes**

The device must be properly configured and have a flushing fluid in the kit. For maximum accuracy of the test results, it is necessary to use the tube system supplied with the device. To connect the pipe system it is recommended to follow the following steps:

- 1) Take the cleaning fluid tubes and waste tubes from the accessories package.
- 2) Connect the tubing system for cleaning fluid and waste according to the color markings on the back of the instrument.
- 3) Connect the fluid level sensors to the cleaning liquid and waste containers.

Picture 2-1. Connection of a system of tubes

Note:

- 1) *The tube system and the power cable must be connected according to the marking on the label.*
- 2) *Do not pull or cut the tubes of the distilled water tank to avoid water ingress under pressure. If necessary, to reduce the pressure, loosen the lid of the container.*

3) Make sure that the distilled water and waste containers are placed on the same surface as the appliance. Avoid bending and straining the tubes, this can lead to both distortion of the results and damage to the device.

2.5 Installation of the keyboard and mouse

If necessary, a keyboard and mouse can be connected to the device. The connection should be performed as follows:

1) Carefully insert the keyboard plug into the corresponding connector on the right side panel of the device, marked with the word "Keyboard".

2) Insert the mouse plug into the corresponding connector on the right side panel of the device, marked with the word "Mouse".

2.6 Installation of the external printer

1) Plug one end of the printer cable into the USB port of the printer.

2) Connect the other end of the printer cable to the USB port of the device.

3) If necessary, install the driver for the connected printer.

4) Connect the printer to a power source using a network cable.

5) Turn on the printer to get started.

2.7 Calibration of the device

Calibration of the device should be carried out using calibration plasmas. There is no need to perform a calibration test daily, but it is mandatory at least once for each new batch of reagent kit that needs calibration. The data obtained as a result of calibration will be used until the data of the new calibration is saved.

Chapter 3 Launch

The list of medical devices and consumables required for use in conjunction with the medical device "Automatic coagulometer "Technology Solution 60" for its intended purpose.

To perform a study of the hemostasis system (testing) in blood plasma according to the principles specified in paragraphs 4.2.1, 4.2.2 and 4.2.3, the device requires reagent kits, calibrators, cuvettes and a washing solution. All components necessary for the operation of the device are specified in section 1.7 of this manual.

Reagent kits.

Reagent kits are necessary for the study of the hemostasis system and are prepared for the operation of the device in accordance with the instructions set out in the relevant instructions for the use of these kits. The procedure for the use of reagents is described in Chapter 6 of this manual.

Disposable cuvettes (see paragraphs 1.2.3 and 1.7 of this manual)

The cuvettes are delivered in packages. K2 cuvettes are installed in tripods for measuring cuvettes, the total number of cuvettes located on board is 72 pcs.

The automatic coagulometer "Technology Solution 60" automatically detects the presence of cuvettes. If no cuvettes were detected in the tripod during operation, the analyzer pauses operation and displays a corresponding message on the screen. In order to be able to share the cuvettes with the coagulometer and the correctness of the results obtained, it is necessary to make sure that there is a special marking on the upper face of each cuvette (see paragraph 1.2.3 of this manual).

Washing solution (see point 1.7 of this manual)

The washing solution is used to clean the sampler during operation. The cleaning cycle is performed according to the factory settings, and additional washing cycles can also be set in the test settings. The procedures for installing and replacing the washing solution are described in Chapter 6 of this manual.

Calibrators

A calibrator is a sample with a known analyte value obtained during production and used for calibration of the test. Depending on the specific method, sets of calibrators can be used (see paragraph 1.7 of this manual), as well as calibrators included in reagent kits. The use of calibrators is performed to obtain the calculated values of the calibrated method. Installation of calibrators on board the analyzer is carried out in the sample compartment using special

cuvettes of the medical device "Disposable cuvette for the line of coagulometers "Technology Solution"", in versions:

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K3 (1000 pcs. in the package).

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K4 (500 pcs. in the package).

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K5 (400 pcs. in the package)

(registration certificate No. RZN 2020/10843 dated June 17, 2020).

Cuvettes K3, K4, K5 are universal for use both when performing calibration procedures and when performing daily in-laboratory quality control. The cuvettes have different capacities: cuvette K3 - 1 ml, cuvette K4 - 1.5 ml, cuvette K5 - 2 ml.

Materials for quality control

Appropriate control materials (plasmas) are used to carry out the quality control procedure. Reagents for quality control are installed in the sample compartment, in special cuvettes. For these purposes, the medical device "Disposable Cuvette for the line of coagulometers "Technology Solution"" is used, in versions:

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K3 (1000 pcs. in the package).

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K4 (500 pcs. in the package).

- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K5 (400 pcs. in the package)

(registration certificate No. RZN 2020/10843 dated June 17, 2020).

Samples

A sample is a portion of the patient's blood plasma. A sample can also be called the material with which the test is performed, in particular calibrators and control plasmas.

The samples are installed in the sample disk. Primary test tubes and sample cups can be used.

3.1 Preparations before launch

Before starting work on the device, the user must perform the following actions:

- 1) Check if there is enough distilled water.
- 2) Empty the waste container.
- 3) Add the required number of measuring cuvettes.
- 4) Empty the tank for the spent cuvettes.
- 5) Make sure that there are no obstacles to the movement of mechanical parts.
- 6) Check whether the device is connected to the network via the UPS.


Note:

When handling liquid waste, always wear disposable medical gloves. After finishing the work, wash your hands in an antiseptic solution.

In contact with blood or other material, infection with pathogenic flora is possible. The issue of disposal of waste materials should be approached seriously, in compliance with all necessary standards.

3.2 Entering the system

After turning on the device and loading the system, a window will appear for logging in:

Picture 3-1. Entering the system

Using the user name “Administrator”, the default password and the password “888888” (you can change the password in the settings after logging in to the system) you can log into

the device as an operator. If you are not using a keyboard, you can click the icon on the monitor, representing the keyboard, to call up the virtual keyboard and enter the password.

3.3 Main menu

After successful login you will see the main menu window:

Picture 3-2. Main menu

The main menu is the main window for all users of the automatic coagulometer "Technology Solution 60", in which users can click the icon to select a specific operation.

When you click on the inscription "Technology Solution 60" in the main menu, a window will appear with information about copyrights and the version of the system.

Picture 3-3. Copyright information and system version

Press "X" to return to the main menu.

3.4 Functional module

All function modules are displayed on the home screen buttons. Their main functions are:

- 1) Measured parameters (Test button): includes the settings of the positions of the reagents, edit the information about the tests and about the urgent sample (details in Chapter 6).
- 2) Quality control parameters (QC button): includes the parameters of materials for quality control, quality control charts, as well as the report and printing of the data obtained after the control (for more details, see Chapter 7)
- 3) Calibration parameters (Calibration button): includes standard settings, standard calibration parameters, manual input of calibration data, as well as a report and printing of data and calibration diagrams (for more details, see chapter 8)
- 4) Parameter Settings (System button): includes information on settings and parameter verification (more details in Chapter 6).
- 5) Query history (Report button): includes querying and printing the history ordered both by patient and by parameter (for more details, see Chapter 9)

- 6) **System Settings (Operator button):** includes information about the hospital, reagents, medical history statistics, combined settings, query log, system parameters, and data loading (for more details, see Chapter 5).

Chapter 4 Settings of parameters

Click the Settings of parameters button in the main menu to select specific options.:

Picture 4-1. Settings of parameters

The first 27 test protocols of the automatic coagulometer "Technology Solution 60" relate to the main built-in parameters. The other three test parameters on this screen are reserved for test protocols created by users as needed. Click Next to see the following test protocols

created by users as needed and select the one you need to open the settings window:

Picture 4-2. Parameter Settings Window

4.1 Settings of positions of parameters

In the settings menu you can set the following test parameters.:

- 1) Test name:

English parameter name - 15 characters.

- 2) Full name:

English parameter name - 30 characters.

- 3) Units of measure:

Selecting the department and viewing information about the medical institution in section 5.4.

Units of measurement of calculation tests:

Some tests require more than one unit to be submitted or require the use of a conversion factor. An example could be a PV test or Prothrombin Time, where to obtain the INR value, the value

MICH. There are several options for obtaining calculation units:

- Ratio, BY (Ratio)
- MHO (INR)

Units of measurement of the calibrated tests:

If the test requires calibration, it is necessary to select the units of measurement in which the calibration curve will be presented. The following units can be selected:

- %,

- g/l (g/l),
- mg/ml (mg/ml),
- mcg/ml (ug/ml),
- mg/dl (mg/dl),
- ng/ml (ng/ml),
- Ed/ml (IU/ml).

The "Calculation Parameter" window includes the possibility of entering reference ranges for each activated parameter. To interpret the test results obtained, it is necessary to refer to the instructions of the set of reagents used on the analyzer (see paragraph 1.7 of this manual).

4) Sensitivity:

The choice of detector sensitivity: low, high and medium.

5) Testing method:

The system includes three methods of testing: coagulation (percentage research method), chromogenic and immunological. Each testing method is divided into 25 types of specific algorithms. For different test methods, you can choose the appropriate test method, for example, for PT as a PT test method, and for FIB - as a FIB test method.

6) Measurement range of test result values. For a fixed percentage of certain parameters, you can set the normal value of the corresponding results and their range; for other parameters, you can take a fixed range of testing time (units of measurement are seconds).

7) Repeat analysis:

Select "Reanalysis". When you select a reanalysis two or more times, the system automatically calculates the average value of the test result; Repeat testing frequency can be selected up to 10 times.

8) Coagulation Point:

Select the stop point of the coagulation measurement for the coagulation method. Choose a stopping point between 2% - 80%, which can vary by an increase of 1%;

Select the start time of the measurement for chromogenic analysis and immunoturbidimetry. Select a time range from 1 second to 600 seconds, which can change by an increase of 1 s..

9) Maximum time:

Choose the maximum clotting time for the coagulation method in the range from 90 to 600 seconds;

Select the measurement stop time for chromogenic analysis and immunoturbidimetry from 1 second to 600 seconds.

10) Sample volume (μL):

Select the volume of sample to be taken in the range of 4 μl ~ 120 μl , which can be changed by an increase of 1 μl . If the diluent is not used, 20 μL or more may be selected.

11) Type of diluent:

Select a diluent and set the reagent parameter. Read more in clause 5.1.

12) Diluent volume:

Select the volume of the diluent to be collected. If the diluents are used only when the standard curve is set, select the diluent name and press 0 μl .

Note:

Ensure that the total volume of the diluent and sample to be withdrawn is in the range of 20 μl ~ 130 μl .

13) Flushing before and after sampling:

Select a flush cycle with a maximum frequency of 9 times before and after sampling, if there is no need for flushing, select "0".

14) S2 sample volume (μL) (second dilution step):

Select the amount of sample to be withdrawn at the second dilution step in the range of 4 μL ~ 120 μL , which can be changed by an increase of 1 μL .

15) Type of diluent (second dilution step):

Select diluent for the second dilution step, then select reagent parameters. Read more in clause 5.1.

16) Diluent volume:

The amount of diluent used in the second dilution step is in the range of 0 μl ~ 200 μl , which can be changed by an increase of 1 μl .

Note: Ensure that the total amount of diluent and sample is in the range of 20 μl ~ 130 μl .

17) Rinsing before and after collection (at the second stage of dilution):

Select a flush cycle with a maximum frequency of 9 times before and after sampling, if there is no need for flushing, select "0".

18) The incubation time of the reagent 1:

Select an interval with a maximum value of 30 seconds between the time the sample starts to warm up and the time it takes to add reagent 1, which can be changed by an increase of 30 seconds. However, if the flush mode is selected immediately after sampling,

each flush cycle will be increased by 30 seconds. In addition, if the flush mode is selected before adding reagent 1, each flush cycle will be increased by another 30 seconds.

19) Name of reagent 1 (R1):

Select the name and parameters of the reagent. Read more in clause 5.1.

20) The amount of reagent taken 1 (R1):

Set the reagent intake volume 1 (R1) between 4 μl ~ 200 μl , which can be changed by increasing 1 μl .

21) Flushing before and after the fence (R1):

Select a flush cycle with a maximum frequency of 9 times before and after collecting reagent 1, if there is no need for flushing, select "0".

22) The incubation time of the reagent 2 (R2):

Select an interval with a maximum value of 30 seconds between the time you add reagent 1 and reagent 2. However, if the flush mode is selected immediately after sampling, each flush cycle will be increased by 30 seconds. In addition, if the flush mode is selected before adding reagent 1, each flush cycle will be increased by another 30 seconds.

23) Name of reagent 2 (R2):

Select the reagent name and reagent parameters. Read more in clause 5.1.

24) The volume of the reagent taken 2 (R2):

Select the amount of reagent to be taken in the range of 4 μl ~ 200 μl , which can be changed by an increase of 1 μl .

25) Flushing before and after collection of reagent (R2):

Select a wash cycle with a maximum frequency of 9 times before and after collecting reagent 2, if there is no need for washing, press "0".

26) The incubation time of the reagent 3:

Select the interval with a maximum value of 30s between the time you add reagent 2 and reagent 3. However, if the flush mode is selected immediately after sampling, each flush cycle will be increased by 30 seconds. In addition, if the flush mode is selected before adding reagent 1, each flush cycle will be increased by another 30 seconds.

27) Name of reagent 3 (R3):

Select the name and parameters of the reagent. Read more in clause 5.1.

28) The volume of the reagent taken 3 (R3):

Select the volume of the reagent to be taken in the range of 4 μl ~ 200 μl , which can be changed by an increase of 1 μl .

29) Washing before and after collecting the reagent (R2):

Select a wash cycle with a maximum frequency of 9 times before and after collecting reagent 2, if there is no need for washing, press "0".

Note: Make sure that the total amount of all solutions in the cuvette is at least 150 µl and not more than 400 µl.

30) Setting up a re-dilution assay:

Samples are retested after a specified multiple dilution. Select the upper and lower limits of the result, as well as the dilution factor. If the test results exceed these values, the sample should be tested again after a specified multiple dilution. After selection, the system automatically starts this function.

31) Configure Reanalysis:

If it is necessary to retest the sample, select the upper and lower limits of the result. If the testing value is outside these limits, the sample will be tested again. After selection, the system automatically starts this function.

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Change measurement results. Formula: $Y = KX + B$; where: X is the measurement result, Y is the variable result. If the changes are not necessary, you can choose the values $K = 1$ and $B = 0$.


33) Click the Save button to save the current parameter value.

34) Click Print. Click the Exit button to exit the current window and return to the main menu.

4.2 Research method parameters

4.2.1 Coagulation method (Percentage study method)

In the optical system of the coagulation method, 90% of the optical scattering + % of the endpoint is used. The light source is an ultra bright red LED with a wavelength of 660 nm, with a high-resolution filter. The light filter directs the light flux of the reaction mixture. The photodiode receives the scattered light flux at an angle of 90° and converts it into an electrical signal. After converting the light flux into a digital signal, the computer processes the data, builds a coagulation reaction curve and determines the coagulation time using the percentage endpoint method.

Picture 4-3. Reaction curve

The clotting time is taken as a percentage.

At the moment when the reagent is just added, the level of scattered light will be 0%, when the coagulation reaction is completed, it will be 100%. Based on the reaction curve, the clotting time can be considered the time when the scattering intensity reaches a certain test percentage. In the figure, the coagulation point is set to 50% and the user can adjust the percentage value of the coagulation point for each time depending on the specific situation. The smaller the value, the shorter the result time will be.

4.2.2 Chromogenic method

The optical system of the chromogenic method uses the optical absorption method. The luminous flux uses an LED, applying a wavelength of 405 nm. The light from the LED passes through the filter and approaches the sample. The photodetector receives light absorption and converts it into an electrical signal. After the analog-to-digital conversion, the computer processes the data, registers changes in optical density per minute and uses the relationship between changes in absorption and concentration (or percentage activity) to check the concentration (or percentage activity) of the test sample.

The chromogenic method is the use of the principle of colorimetry to measure the absorption of light (405 nm) by a liquid in a cuvette. The amount of light that reaches the photodetector is converted into an electrical signal that is proportional to the activity of the enzyme.

Tests that are based on this method: FXIII, AT3, Plg, PCch, Hep, Hep2, a2 AP.

4.2.3 Immunoturbidimetric method

The optical system of the immunoturbidimetry method uses the optical absorption method. The luminous flux uses an LED, applying a wavelength of 575 nm. The light from the LED passes through the filter and approaches the sample. The photodetector receives light absorption and converts it into an electrical signal. After the analog-to-digital conversion, the computer processes the data, registers changes in optical density per minute and uses the relationship between changes in absorption and concentration (or percentage activity) to check the concentration (or percentage activity) of the test sample.

The immunoturbidimetric method is based on measuring the absorption and scattering of monochromatic light by suspended particles formed during the antigen-antibody reaction. There is a change in the optical density of the medium detected photometrically. After plasma incubation for some time, a stabilizing reagent and an antigen-sensitive reagent are added. Then the reaction mixture is placed in a measuring channel with a wavelength (575/660/800 nm, depending on the parameters specified in the test protocol) and the change in light absorption caused by the antigen-antibody reaction is photometrically recorded.

Tests that are based on this method: D-dimer, FDP, PSch, vWF Ag, FM.

Chapter 5 System parameters

In the main menu, click the System Parameters button to see the system settings.

Picture 5-1. System parameters

5.1 Reagent parameters

In the “System Parameters” window, select the Reagent button to enter the reagent parameter settings window.

Picture 5-2. Reagent Settings

In the Reagents window you can select the necessary information about the parameters in the corresponding section. Be sure to complete the following items:

- 1) Name of the reagent
Reagent name in English / Russian, no more than 8 characters.
- 2) Full name
Reagent name in English / Russian, not more than 15 characters.
- 3) Series
Reagent group number, not more than 15 characters.
- 4) Expiration date
Reagent Use Date
- 5) Reagent Type

Type of reagent selected. The system has three types programmed: reagent, washing solution and diluent.

Click the Add button to clear the input field while the cursor is in the reagent name input field.

Click the Save button to save the current reagent settings to memory and the system will automatically save the reagent information according to its type.

Click the Remove button to remove the reagent from the list.


Click Cancel to exit the current window menu and return to the system settings menu.

5.2 Setting up work sheet

In the "Worksheet" window you can select a test group (up to 10 items) to create a worksheet. The system supports three worksheets, each of which can be selected with an independent test group and information on the position of the reagent.

5.2.1 Choice of group of the test

In the "System Parameters" window, select Worksheet Settings to enter the test group settings in the worksheet:


Picture 5-3. Worksheet options

- 1) Select the test group list from the **Select test** list and click the > button to add the test to the worksheet. If such a test has already been added, the operation will not be performed. Press the < button to delete the selected test.
- 2) Click the Save button to save the current worksheet information or the Clear button to delete all the tests.
- 3) Press the Back button to return to the system settings menu..

5.2.2 Group parameters

Click the Group button to enter the Group Settings menu.

Picture 5-4. Group settings

The worksheet groups settings window contains five tests for each of the worksheet groups that will be used to select test parameters of the sample in the group (for more details, see chapter 6.2).

The setup process:

- 1) Select a worksheet.
- 2) Select a group.
- 3) Click on the group name, then you can enter or change the group name.
- 4) Select tests from the list in the **Select test** tab and click the > button to add all selected tests to **List 1**. Click the >> button to add all items in the **Available group options** tab to the **Selected group options** list.
- 5) Select a test from the list of selected group tests, and press the < button to remove selected items from the list, press the << button to clear the list of all available group tests.
- 6) Click the Save Group button to save the selected settings of the current worksheet. Click the Clear Group button to clear the current worksheet settings.
- 7) Click Back to return to the Worksheet Options menu..

Note: If you change the available parameters in the worksheet to which they belong, this will automatically change the parameters of the group.

5.3 The analysis with the autofilter

In the "System Parameters" menu, select the Analysis with AutoFilter button.

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Picture 5-5. AutoFilter Analysis


In this menu, if the user selects an analysis with an autofilter, the system automatically starts performing the corresponding test. For example, if you select the APTT test, if the autofilter value [PT] <0.93, there will be a corresponding result: when the PT test result is <0.93, the system automatically performs the APTT test.

The setup process:

- 1) Click the Filter Parameter button to select the appropriate test in the pop-up window. The autofilter data entry line will show [XXX], where XXX is the selected parameter.
- 2) Click the Operating Parameters button; as a result, the line will display characters for input into the auto filter, such as >, >=, <, <=, OR, AND, and others.
- 3) Repeat steps 1 and 2 to change the filter formula.
- 4) Click the Filter Parameter button to select a test that will automatically appear in a pop-up window.
- 5) Click the Save button to save the current settings or the Delete button to delete the selected entry in the autofilter sheet.

5.4 Information about the hospital

In the "System Parameters" menu, click on the information about the hospital.

Picture 5-6. Settings of the hospital information

Information on the medical institution includes data on the diagnosis, separation, verification of results, etc. (up to 100 entries for each item). To avoid re-entering when editing patient information, you can use information about the patient's illness, department and its results directly from this window.

The setup process:

- 1) Select "Diagnosis", the current parameters will be shown in the list of available parameters on the right.
- 2) Enter the name of the parameter in the Name column.
- 3) Click the Save button to save the current settings.
- 4) Click the Delete button to select and remove a specific item from the list of available options.
- 5) Repeat steps 2-4 to enter the remaining information about the hospital.
- 6) Press the Back button to exit the current menu and return to the System Settings menu.

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
Picture 5-7 System Parameters

- 1) Hospital name: entered by the user, no more than 30 characters.
- 2) Date and time: time and date should be selected according to real time. If you enter incorrect data, the results may be incorrect.
- 3) Printer: The "Technology Solution 60" analyzer model is equipped with a built-in printer whose settings are determined by user preferences.
- 4) Touch Screen: Sensor Setup:
 - a) Select the Save button to save the current settings.
 - b) Click the Touch Screen button to configure the sensor.
 - c) Click the Password button to change the password.

Note: If the printer options do not match the actual configuration, data may not be printed correctly.

5.6 Additional functions

In the "System Parameters" menu, click the Advanced Functions button to see a list of system parameter settings.

Picture 5-8. List of additional features

The system offers two additional functions - monitoring the volume of reagent liquid and analyzing with an auto filter.

5.6.1 Monitoring of volume of liquid of reagent

Allows you to use the monitoring of the volume of liquid reagent during operation.

5.6.2 The analysis with the autofilter

When the autofilter function is present, the system will automatically execute the specified parameters.

Note: Additional features are not available by default.

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Picture 5-9. Query Log


Each system log entry contains information about the user, operation, date of operation, time of operation and various kinds of changes. The operations that the coagulometer “Technology Solution 60” records include startup, shutdown, errors and warnings. These records will be automatically created during the entire time the machine is running. The system log database contains about 1000 records, and the system will automatically display all daily records by default.

The setup process:

- 1) Set the date range of the request.
- 2) Click the Request button to select the date range of the request and the request of the system log.
- 3) Select an entry to delete from the system log list and click the Delete button and select “Delete All” to clear all entries.
- 4) Click the Cancel button to exit the current menu and return to the system settings menu.

5.8 Sending data

In the System Parameters menu, click the Send Data button to enter the corresponding menu.

Picture 5-10. Sending data


Sending data is the transfer of analysis data from the coagulometer “Technology Solution 60” to the laboratory information system for further study.

The setup process:

- 1) 1) Select the date range of the request.
- 2) 2) Click Submit to submit the required test data.
- 3) 3) Click Back to return to the System Parameters menu..

Chapter 6 Test

Click the Test button in the main menu to enter the sample test setup menu.

Picture 6-3. Reaction Block Settings

- A) Date and time: displays the date and time of the current test.
- B) Worksheet: displays the name of the current worksheet.
- C) Temperature: displays the temperature of the test position. Click the Temperature button to open the full system temperature monitoring window. Pressing again closes this window.
- D) System Status: Displays the current status of the system test.
- E) Sample block: there are a total of three racks in the apparatus, each of which can hold 9 samples. Racks have serial numbers: S1-S2-S3.
- E) Sample number in the stand: shows the number of the current sample in the stand. Each such number is unique and is not repeated throughout the day. For example, 002-05 means the number of the fifth sample of the second rack from the beginning of the working day.
- G) Sample number. ID: identification number of the current sample, not more than 20 characters.
- H) Sample Status: shows the status of the current sample.
- I) Type of sample: indicates the type of sample (patient sample and quality control materials)

K) Rack information sheet: displays the information of the current rack compartment.


6.1 The position of the reagent

To carry out testing using the medical device "Automatic Coagulometer "Technology Solution 60" for its intended purpose, it is necessary to use the medical devices specified in paragraph 1.7 of this manual. The Reagent tab is used to place reagents on board the analyzer and identify them in the software. The Reagent tab displays a map of the location of reagents in the reagent tripod. Both reagents and the washing solution, diluents and buffers (paragraph 1.7 of this manual) are placed in the reagent stand, each of which has a different color indication.

The procedure for identifying the reagent on board the analyzer is performed in two ways: manually and automatically (clause 6.4 of this manual).

6.1.1 Setting the position of the reagent

In the "Test" menu, click the Reagent button to enter the reagent position parameters.

Picture 6-2. Reagent Position Settings

To install the required reagent, you need to select the required reagent in the list on the right, then the position on the reagents tripod.

In the reagent position settings menu, you can set the position of the reagent required in the test parameters of each group.

The setup process:

- 1) Select **List 1**, after which all reagents required for the tests in this list will be displayed in the list of required reagents on the right side of the screen.
- 2) Select a reagent from the list of required reagents, and select the position in which it should be placed, after which this position will be highlighted in different colors depending on the type of reagent.
- 3) Click the Delete button and select the reagent item you want to clean. After that, this position will be highlighted in gray, indicating that the operation has been completed successfully. Click the Delete All button to clear all current reagent positions.
- 4) Click the Test button, after which the test parameters of this worksheet will be displayed in the list of test parameters on the right side of the screen.
- 5) Select **List 2 / List 3**, and repeat steps 1-4.
- 6) If a sample cuvette is used, it is necessary to select this type as the cuvette for samples after selecting the reagent position, and then press the Select Capacity button.

Automatic identification of the reagent (barcode scanning)

If there is an external barcode scanner, it is possible to identify the reagent automatically. To do this, enter the "Reagent" menu from the "Analysis" menu. With an external barcode scanner, scan the barcode located on the label of the reagent bottle, and then select the desired position in the menu by pressing. When the function of scanning reagent barcodes to install reagents is enabled, the ability to manually install reagents is blocked.

The procedure for replacing the reagent and the washing solution

Depending on the situation, the reagent can be replaced both during waiting and during sample testing. When the analyzer is in standby mode, it is necessary to get a bottle with a reagent that requires replacement and put a bottle with a new reagent in its place. To detect the reagent level, the analyzer will perform a check before performing the analysis automatically. To replace the reagent, you need to press the "Suspend test" button, then wait for the message that the work is suspended and replace the reagent. Then click the "Start" button and the analyzer will continue working from the moment when it was suspended.

6.1.2 Choose the list.

Select the current list to test a test for a group of tests on this sheet (up to 10 tests).

Final operation:

- 1) Select the current list. Select a worksheet and click Save.

6.2 Reactionary block

The reaction block is designed to set up information about the patient's sample, such as: type of sample, tests performed, type of sample container, and so on.

6.2.1 Information of the reactionary block

In the “Test” menu, click the Reaction Block button to enter the corresponding settings menu.

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A) System Test Status: Displays the status of the current system test (normal testing or urgent sample mode).

B) List: displays the name of the active list.

B) Type of sample: shows information about the patient's sample (green color), urgent sample mode (red color) and quality control materials (yellow).

D) Reaction unit: there are a total of three racks in the apparatus, each of which can hold 9 samples. Serial number: S1-S2-S3.

E) Parameters list: includes a list of all tested parameters of the current list (for more details on the list settings and the selection of the required list and the corresponding reagent position, see clauses 4.2 and 5.1).

E) Combined list: includes all combined test settings in the current list (for more details, see paragraph 4.2). The combined list should contain no more than 10 parameters. It serves to avoid repetition of operations, and also effectively accelerates the work of the reaction unit.

G) Selection of a list of tests: the choice of tests for which a specific sample should be examined.

H) Rack number: shows the number of the current sample holder. Each such number changes daily and is not repeated. For example, 002-05 means the number of the fifth sample of the second rack.

I) Type of sample: indicates the type of sample, including the sample of the individual patient and quality control materials.

K) Sample number. ID: identification number of the current sample, not more than 20 characters.

L) Type tubes for sample storage.

6.2.2 Settings of the reactionary block

The setup process:

- 1) In the Reaction Block tab, select a rack for the sample. This sample will be automatically assigned an identification number displayed on the screen. (The identification number is assigned according to the principle year + month + day + holder number. For example, the number 2018092500203 means September 25, 2018, the third compartment of the second reaction unit).
- 2) Select the Sample Type tab and select the type of this sample rack.
- 3) Users can manually change the sample number. However, if the sample number is already set, it is no longer subject to change.
- 4) On the Test List tab, select the item you want to add. Click the >> Add Item button and the added item will be displayed in the list of selected tests on the right side of the screen (for more details see paragraph 5.2).
- 5) In the Combined List tab, select the group you want to add. Click >> Add group and all test parameters of this group will be displayed in the list of selected tests on the right side of the screen (for more details see paragraph 5.2).
- 6) In the Selected tests list tab, select the item to be deleted, click <Single delete to delete the selected item.
- 7) Click << Delete All to delete the entire list of tests.
- 8) Click Save to save the current settings.
- 9) To enter information for the other reaction chamber, repeat steps 1 through 8.
- 10) Select a rack with a sample, press the Delete button to remove all information from this rack.
- 11) Press the Exit button to exit the current menu and return to the "Test" menu. This menu will display information about all current samples that are in the reaction block, information about which was successfully entered.

6.3 Test of sample

6.3.1 Preparation for a research


- 1) Select the list list and the position of the reagent.
- 2) Enter the sample information.
- 3) Check to see if the position of the reagent matches the settings of the worksheet.
- 4) Ensure a sufficient amount of reagent.
- 5) Make sure there are enough cuvettes for testing.
- 6) Check that the waste container is empty.

- 7) Check the contents of the cleaning tank. If it is empty, please fill.

6.3.2 Barcode scanning (if a scanner is available)

To identify and safely combine medical devices (reagents), it is required to scan the barcode applied to each reagent vial.


In the analysis menu, click the "Tripod" button, in the window that opens, See Figure 6-4, select the tripod that contains the samples intended for identification. Then click the "Start" button.

Picture 6-4. Choosing a tripod for scanning barcode samples

After the "Scan" symbol appears (Fig. 6-5), start sliding the tripod into the device. The scanner will sequentially scan all the installed samples in the tripod. After the scan is completed, all the read sample IDs will be displayed in the list of entered patients.

The success of automatic barcode detection depends on the correct location of the barcode bottle in the reagent stand. The label with the barcode should be visible through the slot of the tripod.

Picture 6-4. Sample barcode scanning procedure

6.3.3 Sample collection

1) Sample collection.

It is necessary to take a venous blood sample. A mixture of venous blood with an anticoagulant in a ratio of 9:1, place a test tube containing 3.8% or 3.2% sodium citrate, then carefully turn it over 10 times to avoid the resulting mixture of clots and bubbles.

2) Disposal of failed samples

All samples containing hemolysis, clots, and other samples unsuitable for further use should be disposed of in accordance with the rules for the disposal of medical waste.

3) Preparation for the test

After taking the sample, it is necessary to separate the plasma as quickly as possible (within 60 minutes), for which the blood is centrifuged at 1200-1600 g for 15 minutes to separate the plasma and shaped elements.

6.3.4 Explanation of the research process

1) The test unit consists of three tripods. (S1-S2-S3). The sample test will be carried out by taking a tripod with samples for an executable module, that is, carried out in the appropriate order S1-S2-S3 (9 samples each).

2) During the testing of the sample of block S1, it will be in the "closed" status, and the other two blocks in the "open" status (the "closed" status means that this tripod with samples has already been sent for testing and users have no opportunity to add additional samples to this tripod; the "open" status says that there is such a possibility).

3) Operations with samples in the block with the status "open".

- Changing sample information: Changing the entered data about the tripod with samples.
- Additional sample: adding a sample to a free position.

- Urgent sample installation: Urgent sample can be sent for testing first. If the placement of samples in the tripod has already been established, the samples for routine testing can be replaced with an urgent sample. If no sample locations have been set, enter the missing data.

Testing procedure

Before starting the testing, it is necessary to perform the preparatory work described in paragraph 3.1 of this manual.

To start the testing procedure, follow the steps listed below:

1. Prepare the primary test tubes by removing the closing caps to avoid collision and possible damage to the sampler.

2. Place the prepared test tubes in the sample compartment. If the test tubes have a barcode, then it is necessary to expand the barcode towards the slot in the sample tripod.

3. In the software interface, specify the necessary testing methods for each sample. If the analyzer is connected to the LIS, then it is necessary to scan barcodes. After successfully scanning the barcodes, press the "LIS" button in the interface to get the task from the LIS database.


4. After all samples have been installed and methods have been obtained or set for them, it is necessary to click the "Start" button located in the lower left corner of the "Analysis" menu.

During the testing process, the analyzer will inform the operator about the status of the samples being performed, as well as display the completed test results. After the analyzer finishes the testing procedure, a message will appear in a separate window on the screen. The procedure for calculating the displayed results is described in Chapter 2 of this manual. To interpret the results obtained, refer to the instructions of the reagent kits (see paragraph 1.7 of this manual).

For the rules for performing repeated measurements, see paragraph 6.3.5. of this manual.

6.3.5 Test Flow Settings


- 1) In the "Test" menu, press the Start button to start the test and enter the initial position of the cuvette.

Picture 6-6. Choosing the starting position of the ditch

2) After selecting the starting position of the cuvette, a window will appear on the screen with the inscription: "Please confirm that the cuvettes / reagents/ samples are placed correctly." After confirmation, you can start testing.


3) The system will take the cuvettes to start the test. Click the Confirm button to exit their current menu and return to the "Test" menu. The testing interface will appear, and the frame with the sample S1 is in the "Locked" state (it turns dark gray).

Picture 6-7. Test status

4) In the "Test" menu, click the Daily Test button in the upper right part of the screen, after which a list of available functions will appear.

5) Click the Add Sample or Add Urgent Sample button to enter the settings of the reaction unit (for more information, see paragraph 6.2).


Picture 6-9. Sample settings during the test.

Note: If the block with samples is in the "closed" status, it is colored gray on the screen, you can save the block settings only in the "open" status.

6) Press the Pause Test button to pause the current test. After pressing this button, all tests will be stopped at the end of the current cycle.

7) In order to view information about the location of samples in the tripod, select a sample entry from the information list of samples in the tripod in the "Test" menu, click the Browse button to display a list of tripods with samples on the screen.

Picture 6-10. Test status

A) Test: shows the test parameter in the test stand.

B) Test Status: shows the status of the test parameter, where R is test readiness, I is incubation, T is the testing process and F is the end of testing.

B) S / OD: Shows test results, time or OD rates.

%, R, INR, DFbg: the result is automatically calculated based on current measurements.

D) Information: displays prompts about the test status for the test parameter.

D) Cancel: Return to the rack information sheet with samples.

E) Curve: Displays the reaction curve of this analysis.

G) Printing: printing the results of a specific sample.

H) Information sheet of the sample compartment. Displays the current rack sample information.

8) Editing basic patient information. Click the Details button to enter the patient master data menu.

Picture 6-11. Editing patient data

9) After testing is completed, you can restart the parameters of any sample compartment. Click the Reaction Block button to enter the corresponding menu and select the test for which you want to retest.

10) Select a rack with a sample and click the Retest button to repeat the selected test.

Picture 6-12. Choosing a test to repeat

Assign a test for retesting, click Retesting to exit the current menu and return to the "Analysis" menu, then click Start to start testing.

Select the item for retesting, press the Retest button to exit the current menu and return to the "Test" menu, then press the Start button to start testing.

Note: Reagents, materials for quality control, as well as the necessary standard materials should be used strictly before the expiration date.

Chapter 7 Quality control


The automatic coagulometer "Technology Solution 60" can display 365 quality control points. With the possibility of simultaneous storage of 12 types of control materials, 27 preset parameters for quality control, as well as 33 programmable parameters for each type of quality control.

In the main menu, click the Quality Control button to select the control parameter:

Picture 7-1. Quality Control Options

Select a specific item to see the selection window of the quality control document.

Picture 7-2. QC document selection window

7.1 Parameters of Materials of Quality control

Select a quality control document, click the Quality Control Materials button to see the required settings and parameters.

Picture 7-3. Quality Control Material Parameters

The parameters of the Quality Control Materials include:

- ◆ The name of the control material
- ◆ Batch number: batch number of control material
- ◆ Quality control level: the Level of quality control materials is divided into Low, Medium and High.
- ◆ Expiration date: means the date by which the materials can be used. If the expiry date of quality control materials has expired, the name of these materials will not be displayed in the quality control information list.

Click the Save button to save the quality control settings.

Click Cancel to return to the window for selecting a quality control document.

7.2 Charts $\bar{X}$

The arithmetic mean value allows you to estimate the real value of the studied parameter of the control sample. The target arithmetic mean of the studied parameter is pre-certified and indicated in the passport of the control plasma for each series of the set.

7.3 Quality Control Charts (Levi Jennings Charts)

Select a quality control document and click the Quality Control Schedule button to enter the corresponding menu

Picture 7-4. Quality Control Graph

- A) Quality Control Test: displays the relevant parameters.
- B) Quality Control Material: shows the name of the control material.
- C) Quality Control Level: The level of quality control materials is divided into low, medium and high.
- D) Party: the series number of the control material.
- D) Expiration Date: Shows the expiration date of the materials.
- E) SD: Standard Deviations
- G) Value: Value
- H) CV: CV level
- K) Quality Control Data: a review of the control data.
- L) Note: notes on the implementation of QC

The KK table shows the calculation of averages, SD (mean square deviation) and CV (%) (coefficient of variation) of Levi-Jennings maps, on which the measurement results of control materials are applied daily or from series to series. Such graphs are constructed for each test and for each level of control material

7.4 Westgard Rules

The quality control settings allow you to set the Westgard rules for tracking the statistical parameters of the results obtained, as well as identifying random or systematic errors. The settings are represented by the following rules:

- R/4s are rejected when 1 control measurement in the group exceeds the average +2s, and the other exceeds the average -2s;

- 4/1s are rejected when 4 consecutive control measurements exceed the same mean +1s or the same mean -1s of the control limit;
- * 10 x deviate when 10 consecutive control measurements fall on one side of the mean.

7.5 Setting quality control parameters

Select the quality control option and click the Quality Control Settings button to enter the corresponding menu:

Picture 7-5 Quality Control Settings

Quality Control options include the following items:

- ◆ Quality Control Test: displays a test whose parameters are monitored.
- ◆ Number of repetitions: means the number of repetitions of the quality control test.
- ◆ Type of Quality Control: By Type of Quality Control can be divided into %, R, INR and dFIB.
- ◆ Internal control: Defines acceptable quality control limits, includes a control limit and multi-standard Quality Control.
- ◆ Press the Exit button to return to the Quality Control Schedule menu.

7.6 Quality Control Limits

The limits of quality control can be divided into control limit and multi-standard quality control.

- 1) For QC limits, it is necessary to set certain parameters. If during the quality control test results exceed the specified values, an alarm will be generated.
 - Upper Limit: The QC data exceeds the upper limit.

- Top indicator: QC data exceeding these norms.
 - Target Value: Quality Control Average
 - Lower score: QC data outside the normal range.
 - Lower limit: QC data outside the lower limit.
- 2) For multi-standard quality control, the following settings must be made.:
- Target Value: The average value for quality control.
 - Standard deviation: The standard deviation of QC.
 - Excess $\pm 2SD$ once: If QC data exceeds $\pm 2SD$, the system will automatically issue a warning that quality control is out of limits.
 - Excess $\pm 3SD$ once: If the QC data exceeds $\pm 3SD$, the system will automatically issue a warning that quality control is out of limits.
 - Excess of $\pm 2SD$ twice: If the QC data exceeds $\pm 2SD$ by two times, the system will automatically issue a warning that quality control is out of limits.
 - Excess $\pm 2SD$ four times: If the QC data exceeds $\pm 2SD$ four times, the system will automatically warn you that quality control is out of limits.
 - Excess $\pm 2SD$ twice in one group: If the QC data exceeds $\pm 2SD$, the system will automatically issue a warning that quality control is out of limits.
 - Excess or underestimate of results by 10 times: if the QC data exceeds or gives too low a result, the system will automatically issue a warning that quality control is out of the set limits.

Conducting quality control:

- Place the control sample on a rack.
 - Register the QC number of the sample.
 - Register the parameters to be analyzed.
 - Press the "Start" button.
- The analyzed data will automatically be saved in the QC-file.


Chapter 8 Calibration

In the main menu, press the Calibrate button to select the test to be calibrated:

Picture 8-1. Calibration Parameters

Select a test to enter the calibration graph menu.


Picture 8-2. Calibration Graphics Menu

- A) Calibrated test: Shows the name of the test being calibrated.
- B) Type of calculated parameter: subdivided into % (concentration), R, INR and dFIB.
- B) Calibration methods: there are a total of four methods: non-linear regression, linear regression, double logarithmic regression, double non-linear logarithmic regression.
- D) Calibration mode: manual (calibration with manual dilution), automatic (calibration with automatic dilution).
- D) Date and time of calibration.
- E) CV%: When recalibrating, the system automatically calculates the CV%.
- G) List of calibration points: there are a total of 6 calibration results with different levels of dilution.
- H) Calibration curve: shows the response curve of each measured parameter.
- I) Recalibration Data List: contains data requiring recalibration.
- K) Information on the reagent: contains all the parameters of the reagent (time, batch number, expiration date, etc.)

8.1 Setting of calibration parameters

In the calibration schedule menu, click the Calibration Settings button to enter the corresponding menu.

Picture 8-3 Calibration Settings

A) Estimated parameter 1: % activity in percent

B) Calculated parameter 2: ratio (PT ratio)

$$\text{PT ratio} = \frac{t}{\text{average value PT}}$$

t = Actual sample clotting time PT,

Average PT = average PT of a normal patient.

C) Estimated parameter 3: INR international normalized ratio

$$\text{MHO} = (\text{PT ratio})^{\text{ISI}}$$

ISI (МНЧ) = International Sensitivity Index, as determined by manufacturers of thromboplastin reagent.

D) Calculated parameter 4: Dfib (PT derivative FIB, calculated fibrinogen).

E) Normal value: normal PT value (patient's normal PT value).

F) Decimal point.

G) Calibration method: there are a total of four methods: non-linear regression, linear regression, double logarithmic regression, double non-linear logarithmic regression.


J) ISI: International Sensitivity Index determined by manufacturers of thromboplastin reagent.

Click the Save button to save the calibration settings.

Click Back to return to the calibration schedule menu.

8.2 Manual entry


In the Calibration Curve menu, click on the “Manual Input” tab to manually enter the data for calibration.

Picture 8-4. Manual entry of calibration data

The setup process:

- 1) In the Calibration Number pop-up window, select the appropriate number of calibration points (up to 6 calibration values).
- 2) Click on the list that appears to see a pop-up window for entering values.

Picture 8-5. Data input

- 3) Enter the number and press the ← button to move the cursor to the left, press the Clear button to delete the entered data, press the Cancel button to exit the current menu.
- 4) Press the OK button to save the changes and exit the current menu.
- 5) To fill in all calibration data, repeat steps 2-4.
- 6) Click on the “Graph” tab to display the curve according to the manually entered data.
- 7) Click the Save Graph button to accept the current curve plot as the last calibration curve.
- 8) Click the Exit button to return to the calibration graph menu.

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Picture 8-6. Calibration Test

- A) Calibrated test: calibrated test.
- B) Type of calculated parameter: divided by% (concentration), R, INR and dFIB.
- B) Calibration methods: there are four calibration methods in total: non-linear regression, linear regression, double logarithmic regression, double non-linear logarithmic regression.
- D) Time and date of calibration.
- E) Concentration of standard material: more details in the reagent manual.

- E) Dilution point number: refers to the number of calibration points.
- G) List of dilution factors for each group: dilution concentration factor.
- H) List of concentrations for each gradient: the system automatically calculates the concentration according to the concentration of the standard material and the corresponding dilution ratio.
- I) A list of test runs: the number of times the test calibration point was repeated.
- K) List of dilution factors: the system has many dilution factors.
- L) Automatic dilution and manual dilution (automatic calibration with pre-calibrated calibrators).
- M) Reagent position: select test list and set reagent position.

8.3.1 Preparation for calibration

Preparation:

- 1) Select the test list and reagent position for which you want to calibrate.
- 2) Check the compatibility of the reagent position and worksheet settings.
- 3) Verify the appropriate amount of reagent.
- 4) Ensure that the reaction cells are suitable for operation.
- 5) Check that the waste container is not filled.
- 6) Check whether it is necessary to fill the tank for cleaning.
- 7) Automatic calibration dilution: the calibration plasma is placed in the first cell of the S1 rack.
- 8) Manual calibration dilution: the prepared calibration plasma will be placed in support S1 in cells one through six.

8.3.2 Automatic calibration dilution

The setup process:

- 1) In the **Concentration Standard Material** window, enter the appropriate concentration.
- 2) In the **List of Dilution Factors** window, select the required coefficient (if necessary, the user can add and remove coefficients on his own).
- 3) Based on the selected dilution factor, the system will automatically apply this factor for each calibration point that will be displayed in the list.
- 4) Based on the ratios of concentration, standard material and dilution factor, the system will automatically calculate the concentration that will be displayed in the concentration list of each calibration point.
- 5) Click the **Test Number** tab to see a pop-up window for entering values. Click OK to save the changes and exit the current menu.
- 6) Repeat step # 5 to enter all stages of the test.
- 7) Calibration plasma is placed in the first cell of rack S1.
- 8) Click the Start button to start testing.
- 9) Upon completion of the test, the calibration graph automatically appears on the screen.
- 10) Click the Exit button to return to the calibration graph menu.

8.3.3 Manual calibration dilution

Setup process

- 1) Select manual dilution.
- 2) Select the number of calibrator standards.
- 3) In the **Concentration** tab of the **Standard Material**, enter the appropriate concentration of the standard material.
- 4) Select the list of test points number to see a pop-up window for entering data, click the OK button to save the changes and exit the current menu.
- 5) Repeat step 5 to enter all test numbers.
- 6) The prepared calibration plasma is placed in the cells of the strut of unit S1 in the first through the sixth.
- 7) Click the Start button to start testing.
- 8) After testing is completed, the calibration graph will automatically display on the screen.
- 9) Click the Exit button to return to the calibration graph menu.

8.4 Printing of calibration data

In the calibration graph menu, press the Print button to print the calibration results.

Chapter 9 Analytical characteristics

9.1 Coefficient of variation (reproducibility).

The coefficient of variation of the method depends on the set of reagents used. For this information, refer to the operating instructions for the appropriate set of reagents, including those specified in paragraph 1.7 of this manual.

9.2 Sensitivity (detection limit).

The sensitivity of the method depends on the set of reagents used. For this information, refer to the operating instructions for the appropriate set of reagents, including those specified in paragraph 1.7 of this manual.

9.3 Linearity (measuring range).

The linearity of the method depends on the set of reagents used. For this information, refer to the operating instructions for the appropriate set of reagents, including those specified in paragraph 1.7 of this manual.

9.4 Influence of interferences.

The measurement process in the presence of various interferences is specified in the operating instructions of the corresponding set of reagents, including those specified in paragraph 1.7 of this manual.

9.5 Limitation of the method.


Restrictions on the use of a medical device are determined by the biological material in which the study is conducted. Human blood plasma stabilized with an anticoagulant (sodium citrate) is used for analysis. Any other material is considered unintended. Also, to clarify the restrictions, refer to the operating instructions for the corresponding set of reagents, including those specified in paragraph 1.7 of this manual.

9.6 Available reference materials.

Reference materials may be control plasmas for the appropriate methods of analysis, including those specified in paragraph 1.7 of this manual.

Chapter 10 History of data


To view the data history, there are two types of query: query by test and query by patient.

Picture 9-1. Request mode

10.1 Inquiry according to the test

Click "By test" and the button "Query" to display a list of results of the query by test.


Picture 9-2. Test Request

The patient information database can contain up to 1000 data points.

10.1.1 Information request according to the test

In the test results menu for the test, click the Request button to see the query window for the test parameter:

Picture 9-3. Request information on the test parameter

- 1) Select the requested parameter in the test parameter pop-up window.
- 2) Set the terms of the request. Select "Data Range" to define the query range. Select Result Range to determine the appropriate range.
- 3) Click the Request button to confirm or the Cancel button to exit the current menu.

10.1.2 Deleting of test results

In the results menu of the query parameters there are the following deletions:

- Delete one record: you need to select one record, and then click Delete to delete the selected item.
- Delete multiple entries: hold down the Ctrl button and select multiple entries at the same time, then click the Delete button.
- Delete all entries: Select "Delete All", then click the Delete button.

10.1.3 Printing of information on test parameter

In the results menu of an information query, the following operations exist for printing information:

- Print one record: select one record and click the Print button to print the selected record.
- Printing multiple entries: hold down the ctrl key and select multiple entries at the same time, then click Print.
- Print All Records: Select "Print All" and click Print.

10.2 Inquiry on the patient

Select By Patient and click the Request button to view the patient information sheet.

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10.2.1 Query of the patient

In the patient information menu, click the Request button to see the patient information request window.

Picture 9-5. Patient request

The setup process:

- 1) Set the query range.
- 2) Select the query terms, such as name, gender, age, separation, diagnosis, sample number, and so on.
- 3) Click the Request button or the Cancel button to exit the current menu..

10.2.2 Change of information of the patient


In the patient information menu, click the Edit button to edit the master data:

Picture 9-6. Entering Patient Data

You can select the patient's gender, separation, diagnosis, and other information from the system database pop-up list. When you finish entering the information, click OK to save or Cancel to cancel the changes.

10.2.3 Viewing of results

After selecting a specific patient, click the View button to see the corresponding window in which the test results for all parameters will be displayed.

Picture 9-7. View Results

In the patient information results menu, click on any column of your choice to sort all the records in ascending or descending order. The screen will display the "▼" or "▲" icons immediately after the name of the current column. If the number of patients is large, sorting may take some time..

10.2.4 Deleting a patient record

- Delete one record: you need to select one record, and then click the Delete button to delete the selected item.
- Delete multiple entries: hold down the Ctrl button and select several entries at the same time, then click the Delete button
- Delete all entries: Select "Delete All", then click the Delete button.

10.2.5 Print the patient's record


In the results menu of an information query, the following operations exist for printing information:

- Print one record: select one record and click the Print button to print the selected record.

- Printing multiple entries: hold down the ctrl button, select multiple entries at the same time, and then click the print button.
- Print All Records: Select “Print All” and click Print.

Chapter 11 Switching off of a system

After testing is completed, the system should be shut down correctly. To do this, press the Shutdown button in the main menu, after which the exit menu appears on the screen.

Picture 10-1. System shutdown

Click the Shut down button to enter the login menu.

Click the Logout button so that the system automatically saves all changes and disconnects. After that, you must turn off the power.

Press the Back button to return to the main menu.

Note: It is recommended to store important information on a CD or any other digital media.

Chapter 12 Device maintenance

12.1 Maintenance of the device performed by the user

Regular maintenance includes daily and monthly manipulations related to the cleaning of sample probes and reagents, mechanical guides, as well as the flushing of the system.

Attention:

All manipulations related to the cleaning of sample probes and reagents, mechanical guides should be carried out only when the device is switched off!

12.1.1 Procedure for daily maintenance

Before starting work on the device, clean the liquid system

Before starting work, perform the "Long-term flushing" procedure. At the same time, it is necessary to make sure that the container for distilled water is filled, open vials with washing liquid are installed in the reagent tray.

After finishing work on the device

• Make sure that the bottle with the washing liquid is open, the distilled water container is filled, then perform the "Long-term washing" procedure.

• Turn off the power of the device.

• Empty the container with the used measuring cuvettes.

• Close the reagent bottles and put them in the refrigerator.

• Empty the waste liquid container.

• Fill the distilled water tank.

• Wipe the sample/reagent probes with an alcohol-free lint-free cloth.

• Wipe the working surfaces of the appliance.

Alcohol characteristics:

Active ingredient: Ethanol;

Volume fraction of ethyl alcohol, % not less than 70

Warning:

(1) The tip of the sample/reagents probe is sharp. Be careful when cleaning to prevent skin damage.

(2) Always wear protective gloves when cleaning the probe

(3) After the cleaning procedure, wash your hands with disinfectant, as there is a risk of infection.

(4) Dispose of waste and liquids in accordance with SanPiN 2.1.3684-21 "Sanitary and epidemiological requirements for the maintenance of urban and rural settlements, water bodies, drinking water and drinking water supply, atmospheric air, soils, residential premises, operation of industrial, public premises, organization and conduct of sanitary and anti-epidemic (preventive) events"

12.1.2 Procedure for monthly maintenance

To clean the distilled water containers and waste containers, rinse the inner walls of the containers with a 20% alcohol solution, then rinse the containers with distilled water.

12.1.3 Performing maintenance when the analyzer is idle

The downtime of the analyzer is considered to be the absence of any of its activities (sample analysis, calibration, etc.) for more than 7 days. After this period, it is necessary to turn on the analyzer and perform the "Long-term flushing" procedure. At the same time, it is necessary to make sure that the distilled water tank is filled, and open vials with washing liquid are installed

in the reagent stand and in the washing zone of the sample probe. In case of further continuation of downtime, it is necessary to carry out the above procedure every 7 days.

Violation of the frequency of work on washing the analyzer when it is idle can lead to the failure of a number of nodes and systems. Hemostasis analyzers are sensitive to prolonged shutdowns and interruptions in operation.

12.2 Maintenance of the device carried out by specialists

This model of the device needs to be serviced every 6 months. Certified specialists who have been trained at the manufacturer's factory and have a valid certificate of admission to this model of the device have access to this type of maintenance

12.3 Replacing the paper of the built-in printer

The built-in printer is located in the lower right part of the device.

- 1) Push the rubber roller away. Pay attention to the temperature so as not to get burned.
- 2) Place a new roll of paper.
- 3) While holding the edge of the paper, press it down with a roller.
- 4) Close the printer cover.

12.4 Adding a washing solution

1) Upon completion of testing, it is necessary to turn the lid of the container with the washing solution counterclockwise.

- 2) Add new liquid to the container.
- 3) Close the lid tightly.


12.5 Replacing the fuse

The fuse is located next to the power switch to the left of the device. There are two fuses for each device.

- 1) Turn off the power.
- 2) Unscrew the fuse clip with a screwdriver counterclockwise.
- 3) Replace with a new fuse type 5A 250 V.
- 4) Place the clamp in the fuse hole.

Warning: Only a fuse of the specified type should be used.

Fuse clip

Picture 11-1. Fuse replacement

12.6 Software update

The automatic coagulometer "Technology Solution 60" uses software version V1.0r, update date 01.2019. If a new version of the software is released, the update is carried out by qualified specialists of the authorized representative in the Russian Federation of Laboratory Solutions LLC.

Chapter 13 System mistakes and their elimination

This Chapter describes common machine system errors and how to troubleshoot them. If you can not fix the problem yourself, please contact the service center of "Laboratory solutions".

Error	Elimination
1) The machine does not start	Check the fuse Check the voltage
2) Error loading program	Turn off the machine and restart after 10 seconds
3) The printer does not start	Check the on/off button position. Check the fuse
4) The printer does not print data	Check the connection.
5) There is no liquid in the tank to clean	Check for leaks capacity to clean Check the pump for correct operation Check the perfusion hose
6) Liquid is poured out of the cleaning container	Check for leaks capacity to clean Check the pump for correct operation Check the perfusion hose
7) Error re-analysis.	Reaction vessel contaminated The volume of the reaction fluid is less than 150 μ l Check whether the reagent is sufficient. Check if the sample is sufficient Poor connection to the line signal Reaction liquid contaminated Expired shelf life of the reagent
8) The volume of liquid taken by the sampling needle is not the same	Check if the needle is installed correctly
9) Quality Control is beyond the established standards.	Check if the test material has expired. Check if you need to adjust the settings Retest using a different method Check the incubation temperature and retest with a new solution
10) Bubbles in the injector	Make sure the piston is closed to prevent leakage. Check the hose connection
11) Contaminated probe	The probe is installed incorrectly, which contributes to leakage Clean the probe or replace with a new one Make sure there is no liquid residue on the probe Make sure the cleaning fluid is not dirty

Install dual flushing in the system settings

12) Incubation chamber
is not heated

Check the temperature settings, set the temperature
parameter to 37 °C.

13) Liquid leakage from
the machine

Check for a hose in the waste container
The waste container is overfilled

2. Symbol description

Sign of the period of use of environmental protection

Meaning of the sign: the Electronic information product contains certain toxic / dangerous substances, the period of use of which is 20 years. You can use it during this period. When this period expires, the product must enter the recycling system.

Chapter 14. Packing, storage, transportation

The device is packed in a plywood box, ensuring the safety of medical products from climatic and mechanical influences during transportation.

Inside the box, the device is surrounded on all sides by foam polymer pads. Also, the strip of foamed polymer must be located inside the analyzer, protecting compartments for reagents and samples, and blocks of the sample taker and node capture ditch. The straps mechanisms move the blocks of the sample taker and the host capturing a ditch must be wrapped with polymer clips to prevent accidental displacement of nodes.

Store the product in dry rooms at a temperature from -20°C to $+55^{\circ}\text{C}$ and relative humidity $\leq 85\%$

Protect from moisture.

Avoid heating or storage in direct sunlight.

The research room must have an operating temperature of $+10^{\circ}\text{C}$ to $+30^{\circ}\text{C}$ and a relative humidity of $\leq 70\%$

If the appliance has cooled down during transport, please wait until it is heated to room temperature before turning it on.

When storing the product, it must be packed in the packaging of the supplier (manufacturer).

It is allowed to transport all types of vehicles in covered vehicles in accordance with the rules of transport, applicable to transport of this type at a temperature of -20°C to $+55^{\circ}\text{C}$ and relative humidity up to 85%.

The device must be transported only in the packaging of the supplier (manufacturer).

It is allowed to transport the product in an additional transport container of the transport company.

The way of packing of transport container with products shall exclude possibility of their movement.

Chapter 15. Manufacturer's guarantees

The service life of the device is 5 years.

An authorized representative in the Russian Federation, Laboratory Solutions LLC, provides a warranty period of operation of the device for 3 years from the date of its commissioning, subject to compliance with the operating rules given in this manual.

The operator must use the device in accordance with this operating manual.

For proper use of the device during commissioning, our representative will conduct training on working on the device.

If there is a malfunction in the device, contact the qualified specialists of our company's representative office in the Russian Federation.

To maintain the warranty on the device, it is necessary to use reagent kits, cuvettes and a washing solution recommended by an authorized representative in the Russian Federation of Laboratory Solutions LLC during the entire warranty period.

Additionally, we draw the attention of the user of the device that the possibility of using medical equipment from one manufacturer together with accessories from another manufacturer is determined by an authorized representative in the Russian Federation of Laboratory Solutions LLC.

In case of violation of the above conditions, the authorized representative in the Russian Federation of Laboratory Solutions LLC disclaims any responsibility for the quality of the results and for the operability of the device, and also has the right to refuse to replace a failed device or any part of it under warranty.

Routine maintenance and maintenance should be carried out during the periods specified in this manual by specialists who have a valid certificate of admission for this model of the device.

The certificate of admission can be issued by an authorized representative in the Russian Federation of Laboratory Solutions LLC, after training on the basis of the head office, service department.

If a malfunction occurs due to improper operation of the device, as well as non-compliance with the above conditions, the warranty is terminated.

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How to use this manual

Thank you for becoming the owner of the automatic coagulometer "Technology Solution 190"! Before using the device, please read the instructions carefully.

To achieve an optimal result, it is necessary to become familiar with the principle of operation of the equipment before embarking on a clinical trial.

This manual will provide information about the installation, daily use, maintenance and other aspects of working with the "Technology Solution 190" automatic coagulometer.

Differences in functions depend on the version of the model configuration. Save all packaging materials - they are needed for storage, transportation and maintenance of the product.

For all technical questions, please contact the distributor.

Warning: *In case if you ignore these precautions when working with the device, this can lead to a distortion of the analysis results, damage to the analyzer's parts and even a threat to the health of the operator.*

Precautionary measures:

Note: *The device is used for clinical studies. A clinical diagnosis, based on the results of the analysis, must be made by a doctor in accordance with the clinical symptoms and a number of other results.*

Notification

Limited Liability Company "Laboratory Solutions" reserves the last right to make changes to the instruction manual.

The illustrations contained in this manual are given as examples, they may not correspond to the product in real performance. Illustrations are not intended to be used for other purposes. Without the written permission of Laboratory Solutions Ltd., the reproduction or modification of this manual is prohibited.

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- assembly, re-adjustment, revision and repair should be carried out by persons who have special permission from Laboratory Solutions LLC;
- the product must be operated in accordance with this manual;
- All electrical devices must comply with international standards.

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Attention

- If the user cannot ensure proper operation / maintenance, the device may fail and even pose a threat to the health of the user.

- Ensure that the instrument is operated under the conditions specified in this manual. Otherwise, malfunction of the device, distortion of the analysis results, damage to the analyzer's parts and injuries are possible.

Packaging symbols

Temperature limit
Transport packaging must be contained within
given temperature conditions.

Fragile, handle carefully. A medical device may
be broken or damaged if not handled with care

Top
Transport packaging must be specified.
side up

Keep away from moisture

Avoid exposure to sunlight.

Do not roll

Symbols on the label:

	Serial number
	Date of manufacture
	Manufacturer
	In vitro diagnostic medical device
	Product compliance with the basic requirements of the directives of the European Union
	The packaging material is made from recycled materials or can be sent for secondary processing.
	Temperature Range

	Danger
	Biological hazard
	Do not touch
	Do not put your hands inside the appliance during operation.

Labeling layouts for the "Technology Solution 190" device

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

Labeling mockups for packaging

Safety Precautions and Safety Tips

For in vitro diagnostic use only. Carefully read the following warning before use and strictly adhere to it.

Attention! Before operating the device, it is necessary to carefully read the safety regulations and strictly observe them.

- If the device emits an unusual smell or smoke, immediately unplug the power supply and unplug the power cord. If the device is used in such a state for a long time, there is a risk of fire, electric shock or injury. Contact your local service representative for an instrument check.
- Take care not to spill blood or reagents on the device and not to drop brackets or paper clips into the device. This may cause a short circuit or smoke. If this happens, unplug the power supply immediately and unplug the power cord. Then contact your local service representative for an instrument check.
- Do not touch the electrical circuits inside the device, especially with wet hands, there is a risk of electric shock.
- Use rubber gloves and the specified tools, parts and components when servicing and inspecting the instrument. To prevent burns caused by engine failure. After completing the job, wash your hands with a disinfectant.
- Be careful when handling samples. Be sure to use latex examination gloves; otherwise infection is possible. If the sample comes into contact with the eyes or a cut on the skin, wash the contact area with plenty of water and consult a doctor immediately.
- Do not touch the internal parts when the device is working, especially moving parts, not adhering to this may result in injury.
- Any micro residue of sample, quality control material, calibration material, etc. on the sampling probe has potential biological hazard. Therefore, do not touch the sampling probe.
- Reagent use and disposal

- Do not allow reagent to enter the skin. In case of contact with skin or clothing, rinse with plenty of clean water to prevent infection.
- If the reagent gets into your eyes, rinse them thoroughly with water and consult a doctor immediately.
- In case of accidental ingestion of the reagent, induce vomiting, drinking plenty of warm saline water, and consult a doctor immediately.
- In case of accidental contact with the reagent, immediately wash the skin thoroughly with water.
- When handling samples or reagents, be sure to use gloves. After completing the job, wash your hands with a disinfectant cleaner to prevent possible infection by pathogenic and other microorganisms.
- When replacing reagent, use precautionary measures and wear gloves.

Power Supply Voltage, Connection, and Grounding

- Ensure that the power source and instrument ground is reliable and stable.
- Do not plug in the power cord to an outlet that does not correspond to a voltage of 100–240 V. Otherwise, a fire or electrical shock may occur.
- When installing, be sure to use the three-wire cable supplied with the device, ensure good grounding and place the device in a place where you can freely turn off the power. Otherwise, fire or electric shock may occur.
- Do not damage the insulation coating of the cable. Do not pull the cable or hang heavy objects with the cable. Otherwise, a short circuit may occur, which may result in electric shock and fire.
- Disconnect power before connecting peripheral equipment. Otherwise, electric shock or damage to the system may occur.

TABLE OF CONTENTS

Packaging symbols	3
Symbols on the label	3
Labeling layouts for the "Technology Solution 190" device	4
Labeling mockups for packaging	6
Safety Precautions and Safety Tips	7
Table of contents	9
Chapter 1. Introduction	13
1.1. Introduction	13
1.1.1. Device name	13
1.1.2. Features	13
1.2. The layout and structure of the device	14
1.2.1. Automatic loading and unloading system ditch	16
1.2.2. Dosing system	17
1.2.3. Thermostatic system	17
1.2.4. Measurement system	17
1.2.5. Cleaning system	18
1.2.6. Data processing system	19
1.2.7. Computer system	19
1.3. Opportunities	19
1.4. Measured Parameters	19
1.5. Technical parameters of the automatic coagulometer "Technology Solution 190"	20
1.6. Description of the composition of the automatic coagulometer "Technology Solution 190"	21
1.7. Description of medical devices intended for use in combination with the medical device "Automatic coagulometer "Technology Solution 190"	23
1.8. PC Configuration	34
Chapter 2. The principle of operation	35
2.1. Coagulation process	35
2.2. Chromogenic method and immunoturbidimetry	36
Chapter 3. Installation	38
3.1. Device unpacking	38
3.1.1. Stages unpacking	38
3.1.2. Conditions of carriage	38
3.2. Installation conditions and operation	39
3.3. Power Supply Requirements	39
3.4. Installation	40
3.5. Calibration	40
Chapter 4. Launch	41
4.1. Preparation before launch	43
4.2. Enabling the system	43
4.3. Indicators	43
4.4. Login	43
4.5. Main menu	44
4.6. Reference	45

Chapter 5. Test Settings	46
5.1. Test parameters	46
5.2. Saving settings and test parameters	48
5.2.1. Calculation of parameters	48
5.2.2. Add test	49
5.2.3. Remove test	50
5.2.4. Print	50
5.2.5. Restore Defaults	50
5.2.6. Test Display List	50
Chapter 6. Reagent Settings	52
6.1. Reagent Parameters	53
6.2. Reagent Information	53
6.3. Preservation	53
6.4. Barcode Scanning (Optional)	53
6.5. Deletion	54
Chapter 7. Sample Settings (Testing)	56
7.1. Sample container	57
7.2. STAT-test	57
7.3. Sample list	57
7.4. Test Information	57
7.5. Sample input	57
7.5.1. Preservation	59
7.5.2. Batch input	59
7.5.3. The number of analyzes	59
7.5.4. Report	60
7.5.5. Print	60
7.5.6. Broadcast	61
7.5.7. Barcode scanning (optional)	62
7.5.8. Skip positions	63
7.5.9. Retest	63
7.5.10. Reaction curve	64
Chapter 8. Quality Control	65
8.1. Quality Control Settings	65
8.2. Quality control materials	65
8.3. Saving	66
8.4. Charts $\bar{X}$	66
8.5. Quality Control Charts (Levi Jennings Charts)	66
8.6. Westgard Rules	67
Chapter 9. Report	68
Chapter 10. Calibration	69
10.1 Calibration Curve	69
10.2 Auto Calibration	69
10.3 Settings	70
10.4 Manual Calibration	71

10.5 Calibration History	72
Chapter 11. Analytical characteristics	73
11.1 Coefficient of variation (reproducibility)	73
11.2 Sensitivity (detection limit)	73
11.3 Linearity (measuring range)	73
11.4 Influence of interferences	73
11.5 Limitation of the method	73
11.6 Available reference materials	73
Chapter 12. System Settings	74
12.1 Additional Functions	74
12.2 Printing	74
12.3 User Defined Parameters	76
12.4 Profile Installation	78
12.4.1 Adding a profile	78
12.4.2 Editing a profile	79
12.4.3 Deletion	80
12.5 Data Dictionary	81
12.5.1 Adding data to the dictionary	81
12.5.3 Deleting a data dictionary	82
12.6 Communication	82
12.6.1 Saving communication settings	82
12.6.2 Connection	83
12.6.3 Disconnection	83
12.7 Barcode Information	83
12.8 User Management	84
12.8.1 Adding a user	85
12.8.2 Editing a user	85
12.8.3 Deleting a user	86
12.8.4 Changing the password	87
12.9 Information about the laboratory	87
12.10 Testing method	88
12.10.1 Saving test method settings	89
12.10.2 Adding a test method	89
12.10.3 Deleting a test method	89
Chapter 13. Service	91
13.1 Maintenance	91
13.1.1 Sleep Mode	91
13.1.2 Exit Sleep Mode	91
13.1.3 Washing the sample probe	91
13.1.4 Washing the reagent probe	91
13.1.5 Bleeding fluid system	92
13.1.6 Purge Fluid System	92
13.1.7 Reloading the sample probe	92
13.1.8 Reloading the reagent probe	92

13.1.9 Reload Capture Block	92
13.1.10 Reload reagent tray	92
13.1.11 Sample Capacity Reload	92
13.1.12 Reloading STAT Capacity Samples	92
13.1.13 Reloading the mechanical part of the device	93
13.1.14 Reboot the device	93
13.2 Data Maintenance	93
13.3 Software Update	94
13.4 Journal	94
13.5 Engineering service	94
Chapter 14. Analyzer Maintenance	96
14.1 Maintenance of the device carried out by the user	96
14.1.1 The order of daily maintenance	96
14.1.2 The procedure for conducting monthly maintenance	97
14.1.3 Performing maintenance when the analyzer is idle	97
14.2 Maintenance of the device carried out by a specialist	97
Chapter 15. System Faults and Their Elimination	98
Chapter 16. Packaging, Storage, Transportation	102
Chapter 17. Manufacturer's Warranty	103

Chapter 1 Introduction

1.1. Introduction

1.1.1. Device name

Automatic coagulometer "Technology Solution 190".

1.1.2. Features

1. The medical device "Automatic Coagulometer "Technology Solution 190" is a fully automatic coagulometer that collects sample data on the principle of a photoelectric sensor based on nephelometry, chromogenic method and immunoturbidimetry. It can be widely used for the diagnosis of hemorrhagic and thrombotic diseases, monitoring and therapeutic monitoring of thrombolytic and anticoagulant therapies, etc.

2. Based on Windows system 7 / 8.1/10/ Ultimate or higher.

3. Fully automatic coagulation analyzer, which allows you to analyze a large number of samples with high performance, with the function of automatic dilution of samples, automatic calibration with the possibility of saving.

4. The possibility of determining fibrinogen by means of a modified Clauss technique.

5. Creation of up to 15 user-configurable profiles, including one or more simultaneously executed tests.

6. Supports continuous addition of STAT samples during analysis.

7. Supports connecting external printers that meet the following criteria:

Classification: matrix, laser, inkjet.

Connection type: USB, wireless.

8. Open system.

9. HL7 medical protocol, which interacts with the laboratory's clinical information system. Supports two-way data exchange with LIS.

10. Arbitrary installation of STAT samples.

11. A large amount of data archive: saving 150,000 calibration curves, 150,000 reports on the results of the study.

12. Generates a variety of reports for printing analysis results.

13. The standard configuration type of the device is desktop. No additional stand is required.

14. Availability of sample scanners and reagents.

The product includes:

1. Automatic coagulometer "Technology Solution 190" - 1 pc.
2. Operation manual - 1 pc.
3. Instrument passport - 1 pc.
4. Capacity for distilled water - 1 pc.
5. Container for discharging liquid waste - 1 pc.
6. Power cord - 1 pc.
7. Wire for ground connection - 1 pc.
8. Fuse - 2 pcs.
9. Hose for waste disposal - 1 pc.

10. LAN cable - 1 pc.
11. Cup for samples 1.5 ml (if necessary) - no more than 100 pcs.
12. Handles for carrying the device (if necessary) - no more than 4 pcs.
13. Micro adapter for samples (if necessary) - no more than 2 pcs.
14. Adapter for sample cups (if necessary) - no more than 2 pcs.
15. Guide plate for used cuvette (if necessary) - no more than 1 pc.
16. Bottle for reagents $\varnothing 23$ mm (if necessary) - no more than 2 pcs.
17. Bottle for reagents $\varnothing 35$ mm (if necessary) - no more than 5 pcs.
18. Bottle for reagents $\varnothing 28$ mm (if necessary) - no more than 2 pcs.
19. Magnetic stirrer (if necessary) - no more than 2 pcs.
20. Label for reagents (if necessary) - no more than 2 pcs.
21. CD with software (if necessary) - no more than 1 pc.
22. Reader for reagents (if necessary) - no more than 1 pc.
23. Reader for samples (if necessary) - no more than 1 pc.

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1.2 The layout and structure of the device



Figure 1-1. The main view of the device.

The layout includes the instrument and PC, the device is operated using a PC. WLAN interface is used for communication between instrument and PC. The device can be mounted on a stand (not included in the scope of delivery). The main parts of the device: a hopper for loading a reaction cell, a sampling system and reagents, a thermostat, a measurement system, a washing system, a data processing system, and a computer system.


View of the front panel "Technology Solution 190"».

Legend: 1 - Front Panel. 2 - Status screen. 3 - Protective cover. 4 - Cover the hopper loading cuvette. 5 - Top panel. 6 - Right sidebar. 7 - On / Off button. 8 - Capacity for waste cuvette. 9 - Secondary switch.


View of the rear panel "Technology Solution 190"».

Legend: 1 - Rear Panel. 2 - Power supply connection point (as well as two fuses). 3 - Grounding. 4 - WLAN interface. 5 - Left sidebar. 6 - Outlet of waste liquid. 7 - Distilled water inlet.

Legend: 1 - Disk samples. 2 - Positions for STAT samples. 3 - The position of cleaning the sampler. 4 - Automatic ditch loading unit. 5 - Sampler. 6 - Reagent probe. 7 - Module measuring system. 8 - Place of discharge of spent cuvettes. 9 - Position cleaning reagent probe. 10 - Frame. 11 - Reagent Disk.

Legend: 1 - Power Module. 2 - Board Module. 3 - Dosing piston sampler and probe reagents. 4 - Liquid valve. 5 - Feed pump. 6 - Valve block waste fluid.

1.2.1 Automatic loading and unloading system ditch

The automatic loading mechanism and the receiving mechanism of the measuring cuvettes carry out the distribution of the cuvette in the measuring unit, taking them into account and discarding the already processed cuvettes into the tank for the used cuvettes.

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1. Determination of the liquid level. If there is not enough reagent or sample, the device will issue a corresponding message and thereby notify the user.
2. Preheating of the reagent. Before adding the reagent, it is heated to 37 ° C in 3-5 seconds.
3. Prevention of horizontal and vertical collision.

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1. Cooling: the temperature of all reagent positions is $\leq 16^{\circ}\text{C}$.
2. Reagent preheating: the temperature of the reagent probe is maintained at 37 ° C. The reagent is preheated in 3-5 seconds.
3. The incubation unit has 18 incubation cells with a constant temperature of 37 ° C.
4. The measuring unit has 13 measuring cells with a constant temperature of 37 ° C.

1.2.4 Measurement system

The measurement system turns on the light flux, collects signals and preprocesses them during the reaction. The measuring system of the device includes three measurement methods: coagulation, chromogenic and immunoturbidimetric methods.

1. Coagulation method: 8 measuring channels are involved: 1, 2, 3, 4, 5, 6, 7 and 8 marked in black, as shown in the figure below. The wavelength is 660 nm.

2. Chromogenic and immunoturbidimetric methods: 5 measuring channels are involved: 1, 2, 3, 4 and 5, shown in the figure below. The following wavelengths are supported: 405 nm, 575 nm, 660 nm, 800 nm.

3. Incubation of the reaction mixture takes place in positions 1-18, marked in red.
4. Positions 1 and 2 are the discharge positions of the spent cuvettes.


Figure 1-2. Top view of the measuring positions of the measurement unit.

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- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 µl) – 1000 pcs. in the package.
- Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 µl) – 1500 pcs. in the package.

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- Volume 700 µl;
- The total height of the cuvette is 30.06 ± 0.15 mm;
- Diameter of the landing part of the cuvette - 10.00 ± 0.1 mm;
- External diameter of the bottom of the cuvette - 7.70 ± 0.1 mm;
- Distance of the measurement area from the mating plane - 21.40 ± 0.1 mm;
- External diameter of the measuring part of the cuvette - 7.75 ± 0.1 mm;
- Wall thickness in the measuring part of the cuvette is not more than - 1.08 ± 0.03 mm;
- The maximum deviation of the cell wall thickness in the measuring part is 0 ... + 0.01 mm.
- Displacement radius – 1.60 ± 0.03 mm.

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1.2.5 Cleaning system

The cleaning system is designed for flushing the liquid line and cleaning the reagent probe and the sampler.

The cleaning system consists of the following components:

1. Washing bath: designed for jet cleaning of the inner and outer walls of the reagent probe and sampler.
2. Distilled water tank: Used to store distilled or deionized water.
3. Waste liquid container: used to collect waste liquid.
4. As a washing and disinfecting reagent for cleaning the dosing systems of the coagulometer, a medical device must be used in accordance with registration certificate No. RZN 2019/9048 dated May 04, 2021 "Washing solution according to technical specifications 21.20.23-001-34065548-2019".

1.2.6 Data processing system

The system processes and analyzes data and performs physical, chemical and biological analysis.

1.2.7 Computer system

The system coordinates various activities and monitors the operation of the device, including the interaction of man and equipment, displays information and performs printing.

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1.4. Measured Parameters

All measured parameters are shown in Table 1-1

Table 1-1

Name	Measured Parameter	Analysis Method	Units
Prothrombin time	PT	Coagulation	s (second)
Activated partial thromboplastin time	APTT	Coagulation	s (second)
Fibrinogen (standard Clauss method)	FIB	Coagulation	g / l
Fibrinogen (modified Clauss method)	FIB MMC	Coagulation	g / l
Thrombin time	TT	Coagulation	s (second)
Lupus Anticoagulant LA1	LA1	Coagulation	s (second)
Lupus Anticoagulant LA2	LA2	Coagulation	s (second)
Coagulation factor	II, V, VII, VIII, IX, X, XI, XII	Coagulation	%
Coagulation factor	XIII	Chromogenic method	%
Antithrombin III	AT3	Chromogenic method	%
Plasminogen	Plg	Chromogenic method	%

Protein C	PCch	Chromogenic method	%
Protein C	PCco	Coagulation	%
Protein S	PSch	Chromogenic method	%
Protein S	PSco	Coagulation	%
Heparin	Hep	Chromogenic method	IU / ml
DOAC	Hep2	Chromogenic method	ng / ml
α 2 antiplasmin	α 2 AP	Chromogenic method	%
Fibrinogen degradation products	FDP	Immunoturbidimetry	mkg / ml
D-Dimer	D-dimer	Immunoturbidimetry	ng / ml
Antigen factor Von Willebrand	vWF Ag	Immunoturbidimetry	%
Fibrin-monomer	FM	Immunoturbidimetry	mkg / ml

1.5. Technical parameters of the automatic coagulometer "Technology Solution 190"

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Testing method	Coagulation process (optical scattering method), chromogenic method and immunoturbidimetry (optical transmission method)
Wavelengths	Coagulation method: 660 nm Chromogenic and immunoturbidimetric methods: 405 nm, 575 nm, 660 nm, and 800 nm.
Measuring Positions	13 positions
Incubation positions	18 positions
Discharge positions for ditches:	2 positions
Performance:	190 tests / hour (PT test)
Sample positions:	50 positions for samples, 5 positions for STAT samples
Reagent positions:	24 cooled positions, 5 non-cooled positions
Maximum load of measuring cuvettes:	300 measuring cuvettes
Nutrition:	a.c.100-240V, 50Hz / 60Hz
Power:	500 watts
Terms of Use:	10 °C ~ 30 °C, relative humidity $\leq$ 70%
Storage conditions:	- 20 °C 55 °C, relative humidity $\leq$ 85%, Pressure: 50 kPa ~ 106 kPa, in a well ventilated room
Interface:	WLAN, RS-232 port, USB port
Safety class for software acc. To EN 62304	A
electrical safety according to EN 61010-1	- Overvoltage category I. - Degree of contamination 1. - Protection class I.
Group of dependence on perceived mechanical influences according to	2 - wearable, portable and mobile, not intended for use when carrying and moving within a stationary room

GOST R 50444-92	
The degree of protection of the analyzer from dangerous penetration of water or solid particles - IP	IP 00

1.6. Description of the composition of the automatic coagulometer "Technology Solution 190"

№ п/п	Description	Assignment	Technical specifications
	Capacity for a distilled water	It is intended for a distilled water	Volume: 20 l. Overall dimensions: Height: 360 mm $\pm$ 1 mm; Width: 280 mm $\pm$ 1 mm; Depth: 280 mm $\pm$ 1 mm.
	Capacity for draining of liquid waste	It is intended for collecting the fulfilled liquid	Volume: 20 l Overall dimensions: Height: 360 mm $\pm$ 1 mm; Width: 280 mm $\pm$ 1 mm; Depth: 280 mm $\pm$ 1 mm.
	Cup for samples of 1.5 ml (if necessary)	Are used for samples or scarce plasmas during test by the device	Diameter: 15 mm $\pm$ 1 mm. Height: 25 mm $\pm$ 1 mm.
	Network cord	Provides electricity feed-in to the device	Length: 1800 $\pm$ 10 mm. Voltage: up to 500 V. Diameter: 7 mm $\pm$ 0.5 mm.
	Wire for connection of grounding	Provides protective grounding of the device	Voltage: up to 600 V Length: 5200 mm $\pm$ 10 mm. Diameter: 4 mm $\pm$ 0.5 mm.
	Handles for carrying the device (if necessary)	Designed to carry the device	Diameter: 28 mm.
	The micro adapter for samples (if necessary)	Designed for installing cups in a tripod for reagents	Diameter: 24 mm $\pm$ 1 mm. Height: 29 mm $\pm$ 1 mm.
	The adapter for cups of samples (if necessary)	It is intended for installation of cups with scarce plasmas in a support for reagents	Diameter: 24 mm $\pm$ 1 mm. Height: 29 mm $\pm$ 1 mm.
	The directing plate for fulfilled a ditch (if necessary)	It is intended for the direction of reset fulfilled a ditch in capacity outside the device	Width: 118 mm.
	Hose for draining of waste	It is intended for replacement of a hose for draining of waste	Length: 500 mm $\pm$ 5 mm. Diameter: 12 mm $\pm$ 0.5 mm.
	Safety lock	Performs protective function, protects an electrical circuit and its elements from overheating and ignition at course of current of high force.	Diameter x Length: 5x20 mm. Current strength: 10 A. Voltage: 250 V.

	LAN cable	It is intended for connection with the managing computer	Length: 2500 mm $\pm$ 10 mm. Diameter: 0.7mm $\pm$ 0.5mm.
	Bottle for mm Ø23 reagents (if necessary)	It is intended for reagent storage	Volume: 25 ml. Overall dimensions: Diameter: 23mm $\pm$ 1mm; Height: 45 mm $\pm$ 1 mm.
	Bottle for mm Ø35 reagents (if necessary)	It is intended for reagent storage	Volume: 30 ml. Overall dimensions: Diameter: 35mm $\pm$ 1mm; Height: 60 mm $\pm$ 1 mm.
	Bottle for mm Ø28 reagents (if necessary)	It is intended for reagent storage	Volume: 20 ml. Overall dimensions: Diameter: 28mm $\pm$ 1mm; Height: 50 mm $\pm$ 1 mm.
	Magnetic mixer (if necessary)	It is intended for hashing of reagents	Length: 15 mm.
	Sticker for reagents (if necessary)	It is used for designation of installation sites of reagents in a support	Length: 13 mm $\pm$ 0.5mm. Width: 5 mm $\pm$ 0.5 mm.
	CD with the software (if necessary)	It is intended for installation of the software managing the analyzer on the computer	Capacity: 700 MB. CD type: Laser.
	The reader for reagents (if necessary)	It is intended for reading of a barcode of reagents	Specifications: • Wavelength: 650 nm. • Power: up to 1 MW. • Laser class: 2. • Maximum distance to the barcode: 90 cm. • Reading speed: 104 $\pm$ 12 barcodes per second. • Power consumption: 95 MW. • System connection: hardware.
	The reader for samples (if necessary)	It is intended for reading of a barcode of samples	Specifications: • Wavelength: 650 nm. • Power: up to 1 MW. • Laser class: 2. • Maximum distance to the barcode: 90 cm. • Reading speed: 104 $\pm$ 12 barcodes per second. • Power consumption: 95 MW. • System connection: hardware.

1.7. Description of medical devices intended for use in combination with the medical device "Automatic coagulometer "Technology Solution 190", not included in the delivery package.

Name	Number and date of issue of the registration certificate	Appointment
A set of reagents for the determination of activated partial thromboplastin time for the line of coagulometers "Technology Solution" (TS-APTV) according to TU 21.20.23-082-42349142-2020	RZN 2020/12933 of December 14, 2020	The medical device "TS-APTV" is designed to perform a basic technique for the study of the hemostasis system - the determination of activated partial (partial) thromboplastin time (APTV / APTT) on the line of semi-automatic and automatic coagulometers "Technology Solution". APTT is determined to evaluate the internal pathway of coagulation activation in blood plasma.
A set of reagents for determining the prothrombin time for the line of coagulometers "Technology Solution" (TS-Thromboplastin) according to TU 21.20.23-083-42349142-2020	RZN 2020/12930 of December 15, 2020	The medical device "TS-Thromboplastin" is designed to determine the prothrombin clotting time (PV) in venous blood plasma, calculate the prothrombin ratio (PO), the prothrombin index (PTI), the international normalized ratio (INR) and the prothrombin index by Quick on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of prothrombin time is carried out when monitoring the treatment with anticoagulants of non-direct action, as well as to detect dysfunction of the external coagulation pathway.
A set of reagents for determining thrombin time for the line of coagulometers "Technology Solution" (TS-Thrombin) on the	RZN 2020/12659 of November 24, 2020	The medical device "TS-Thrombin" is designed to determine the thrombin time on the line of semi-automatic and

TU 21.20.23-084-42349142-2020		automatic coagulometers "Technology Solution". Determination of thrombin time is carried out to assess the final stage of blood plasma coagulation during prophylactic or therapeutic administration of anticoagulants, in connection with arterial or venous thrombosis (myocardial infarction, deep vein thrombosis, pulmonary embolism, etc.), as well as in the detection of congenital or acquired (cirrhosis of the liver) dysfibrinogenemia.
A set of reagents for determining the concentration of fibrinogen for the line of coagulometers "Technology Solution" (TS-Fibrinogen) on the TU 21.20.23-085-42349142-2020	RU No. RZN 2020/12640 of November 30, 2020	The medical device "TS-Fibrinogen" is intended for quantitative determination of fibrinogen in blood plasma without preliminary dilution of the studied plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". Determination of fibrinogen concentration is carried out with bleeding of various etiologies (during pregnancy and childbirth, against the background of fibrinolytic therapy, etc.), with inflammatory, autoimmune, infectious and septic diseases (pyelonephritis, pneumonia, peritonitis, sepsis, systemic lupus erythematosus, etc.).
A set of calibrators for determining the concentration of fibrinogen for the line of coagulometers "Technology Solution" (TS-Fibrinogen calibrator) according to TU 21.20.23-086-42349142-2020	RZN 2020/12691 of November 30, 2020	The medical device "TS-Fibrinogen-calibrator" is designed to obtain calibration values of clotting time when determining the concentration of fibrinogen in blood plasma by the modified Clauss method without preliminary dilution of the plasma under study on the line of semi-automatic and automatic coagulometers "Technology Solution".

		Solution". The medical device "TS-Fibrinogen-calibrator" is intended only for the set "TS-Fibrinogen".
A set of reagents for the determination of antithrombin for the line of coagulometers "Technology Solution" (TS-Antithrombin) on the TU 21.20.23-087-42349142-2020	RZN 2020/12686 of November 30, 2020	The medical device "TS-Antithrombin" is designed to determine the activity (as a percentage of the norm) of the physiological anticoagulant – antithrombin (AT) on the line of coagulometers "Technology Solution". Determination of AT is carried out to diagnose congenital and acquired deficiency of this anticoagulant (DIC syndrome, various types of shock, sepsis, obstetric pathology, severe trauma, etc.), as well as with prolonged or high-dose use of heparins.
A set of reagents for the quantitative determination of D-dimer in blood plasma for the line of coagulometers "Technology Solution" (TS-D-dimer) on the TU 21.20.23-088-42349142-2020	RZN 2020/12270 of November 14, 2020	The medical device "TS-D-dimer" is intended for the quantitative determination of D-dimer in human blood plasma by the immunoturbidimetric method at a wavelength of 570-800 nm on the line of coagulometers "Technology Solution". D-dimer is one of the most important markers of activation of blood clotting and fibrinolysis, since it is formed in the bloodstream as a result of the formation and subsequent lysis of stabilized fibrin by plasmin. A high level of D-dimer will be registered in patients with thrombosis, thromboembolism, DIC syndrome of various genesis.
TS-D-dimer control	RZN 2020/12571 of November 17, 2020	The medical device "TS-D-dimer control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of

		D-dimer in blood plasma for the line of coagulometers "Technology Solution" (TS-D-dimer)" on the line of coagulometers "Technology Solution".
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS-calibrator) on the TU 21.20.23-090-42349142-2020	RZN 2020/12354 of October 28, 2020	The medical device "TS-calibrator" is used to obtain calibration values and build calibration curves for the study of the hemostasis system on the line of semi-automatic and automatic coagulometers "Technology Solution". The use of calibration plasma frees in most cases from the need to receive fresh normal blood plasma from healthy people. Calibration plasma "TS-calibrator" is certified for constructing calibration curves and obtaining calibration values when determining the following indicators: - APTV/ACTV; - prothrombin time; - prothrombin index by Quick; - thrombin time; - antithrombin; - plasminogen; - protein C; - Willebrand factor antigen level (vWF:Ag); - the level of protein S (free) antigen.
A set of control plasmas for the line of coagulometers "Technology Solution" (TS-control) according to TU 21.20.23-091-42349142-2020	RZN 2020/12097 of September 30, 2020	The medical device "TS-control" is designed for in-laboratory quality control (Quality Control) in the study of the hemostasis system on the line of semi-automatic and automatic coagulometers "Technology Solution". The set is certified in the normal

		and pathological ranges according to the following indicators: - APTV/ACTV; - prothrombin time; - international normalized attitude (INR); - prothrombin index by Quick; - thrombin time; - fibrinogen (modified by Klaus method); - antithrombin; - plasminogen; - protein C; - factor VIII; - factor IX; - Willebrand factor antigen level (vWF:Ag); - the level of protein S (free) antigen.
A set of reagents for determining the activity of factor VIII for the line of coagulometers "Technology Solution" (TS-Factor VIII) according to TU 21.20.23-092-42349142-2020	RZN 2020/13132 of January 12, 2021	The medical device "TS-Factor VIII" is designed to determine the activity of factor VIII in blood plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of factor VIII is used for the diagnosis of hemophilia A, for monitoring the replacement therapy of patients with hemophilia A with factor VIII compounds, as well as for the diagnosis of thrombophilia caused by an increased level of factor VIII.
A set of reagents for determining the activity of factor IX for the line of coagulometers "Technology Solution" (TS-Factor IX) according to TU 21.20.23-093-42349142-2020	RZN 2021/13310 of January 28, 2021	The medical device "TS-Factor IX" is designed to determine the activity of factor IX in blood plasma on the line of semi-automatic and automatic coagulometers "Technology Solution". The determination of factor IX is used to diagnose hemophilia B, to control the

		replacement therapy of patients with hemophilia B with factor IX concentrates.
A set of reagents for the determination of plasminogen for the line of coagulometers "Technology Solution" (TS-Plasminogen) on the TU 21.20.23-095-42349142-2020	RZN 2021/13401 of February 08, 2021	The medical device "TS-Plasminogen" is designed to determine the activity (as a percentage of the norm) of the main component of the fibrinolytic system – plasminogen on the line of coagulometers "Technology Solution". Determination of plasminogen is used in the diagnosis of DIC syndrome and thrombophilia, control of treatment with drugs that activate fibrinolysis.
A set of reagents for the determination of protein C for the line of coagulometers "Technology Solution" (TS-Protein C) according to TU 21.20.23-096-42349142-2020	RZN 2021/13611 of March 02, 2021	The medical device "TS-Protein C" is designed to determine the activity (as a percentage of the norm) of protein C in blood plasma on the line of coagulometers "Technology Solution". Protein C is one of the most important physiological anticoagulants synthesized in the liver with the participation of vitamin K. In its active form, it destroys coagulation coagulation factors VIIIa and Va. Hereditary protein C deficiency leads to thrombosis and thromboembolism. Newborns with protein C deficiency may develop fulminant purpura.
A set of reagents for the determination of lupus anticoagulant for the line of coagulometers "Technology Solution" (TS-Lupus) according to TU 21.20.23-094-42349142-2020	RZN 2021/14378 of May 20, 2021	The medical device "TS-Lupus" is intended for the determination of lupus anticoagulant (VA) on the line of semi-automatic and automatic coagulometers "Technology Solution" in sodium citrate-stabilized blood plasma. In blood plasma, VA binds to negatively charged phospholipids and protein-phospholipid

		complexes and inhibits the activation and interaction of plasma coagulation factors. These disorders are most clearly detected in phospholipid-dependent coagulation tests. The presence of VA in plasma is accompanied by recurrent thrombosis of veins and arteries, disorders of cerebral circulation (headaches, fainting, dynamic disorders of cerebral circulation, paresis, episynndrome, visual disturbances, etc.), fetoplacental insufficiency, habitual miscarriage (miscarriages, fetal death of the fetus), thrombocytopenia, less often – microcirculatory bleeding, polyallergia, others immune disorders, a tendency to develop DIC syndrome.
A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa) on the TU 21.20.23-102-42349142-2021	RZN 2021/15119 of August 18, 2021	The medical device "TS-Anti-Xa" is intended for quantitative determination of the anti-Xa activity of direct anticoagulants in blood plasma by chromogenic method on the line of coagulometers "Technology Solution" when monitoring the therapy with anticoagulants, such as unfractionated and low molecular weight heparins, direct oral anticoagulants (DOAC) (rivaroxaban, apixaban, fondaparinux and others).
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS-Heparin calibrator) on the TU 21.20.23-103-42349142-2021	RZN 2021/15136 of August 20, 2021	The medical device "TS-Heparin calibrator" is designed to obtain calibration values for determining the anti-Xa activity of direct anticoagulants (unfractionated and low molecular weight heparins) in blood plasma by chromogenic method on the line of

		coagulometers "Technology Solution". The medical device "TS-Heparin calibrator" is intended only for the "TS-Anti-Xa" kit.
A set of control plasmas for the line of coagulometers "Technology Solution" (TS-Heparin control) according to TU 21.20.23-105-42349142-2021	RZN 2021/15104 of August 16, 2021	The medical device "TS-Heparin-control" is intended for in-laboratory quality control of a medical device "A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa)" when determining the anti-Xa activity of direct anticoagulants (unfractionated heparin (NG) and low molecular weight heparins (NMH)) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution".
A set of calibration plasmas for the line of coagulometers "Technology Solution" (TS- DOAC -calibrator) on the OTHER 21.20.23-104-42349142-2021	RZN 2021/15411 of September 23, 2021	The medical device "TS- DOAC -calibrator" is designed to obtain calibration values for determining the anti-Xa activity of direct oral anticoagulants (DOAC) (apixaban, rivaroxaban, fondaparinux and others) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution". The medical device "TS- DOAC -calibrator" is intended only for the set "TS-Anti-Xa".
A set of control plasmas for the line of coagulometers "Technology Solution" (TS- DOAC -control) according to TU 21.20.23-106-42349142-2021	RZN 2021/15147 of August 20, 2021	The medical device "TS- DOAC -control" is intended for in-laboratory quality control of a medical device "A set of reagents for determining the anti-Xa activity of direct anticoagulants for the line of coagulometers "Technology Solution" (TS-Anti-Xa)" when determining

		the anti-Xa activity of direct oral anticoagulants (DOAC) (apixaban, rivaroxaban, fondaparinux and others) in blood plasma by chromogenic method on the line of coagulometers "Technology Solution".
A set of reagents for the determination of the Willebrand factor antigen for the line of automatic coagulometers "Technology Solution" (TS-Willebrand Ag) according to TU 21.20.23-107-42349142-2021	RZN 2021/15490 of October 04, 2021	The medical device "TS-Willebrand Ag" is intended for the quantitative determination of the Willebrand factor antigen (vWF:Ag) in blood plasma. Willebrand factor (vWF) is a large multimeric glycoprotein whose main functions are participation in platelet adhesion by binding it to glycoprotein Ib, as well as protection of antihemophilic globulin from proteolysis. Definition of vWF:Ag is one of the most popular methods of diagnosing hemorrhagic disease, which is based on a deficiency or functional inferiority of vWF (Willebrand's disease, Willebrand syndrome). In addition, it is advisable to use this method to monitor the effectiveness of the treatment of this pathology with Willebrand factor concentrates.
A set of reagents for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer) on the TU 21.20.23-113-42349142-2021	RZN 2022/16545 of February 15, 2022	The medical device "TS-Fibrin-monomer" is intended for the quantitative determination of fibrin monomers (FM) in human blood plasma either in free form or as part of soluble fibrin-monomer complexes (RFMC) by immunoturbidimetry on coagulometers "Technology Solution". The definition of FM is important in various variants of intravascular blood clotting (thrombosis of the main vessels, PE, DIC syndrome), including the

		use of anticoagulants.
A set of control plasmas for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer control) according to TU 21.20.23-115-42349142-2021	RZN 2022/19288 of December 30, 2022	The medical device "TS-Fibrin-monomer control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer)". The medical device "TS-Fibrin-monomer control" is intended only for the set "TS-Fibrin-monomer".
A set of calibration plasmas for quantitative determination of FM on coagulometers "Technology Solution" (TS-Fibrin-monomer calibrator) according to TU 21.20.23-117-42349142-2021	RZN 2022/19253 of December 27, 2022	The medical device "TS-Fibrin-monomer calibrator" is designed to obtain calibration values for the quantitative determination of fibrin-monomer (FM) with a set of "TS-Fibrin-monomer" on coagulometers "Technology Solution".
A set of reagents for quantitative determination of FDP on coagulometers "Technology Solution" (TS- FDP) according to TU 21.20.23-114-42349142-2021	RZN 2022/19267 of December 28, 2022	The medical device "TS- FDP " is designed for the quantitative determination of fibrinogen and fibrin degradation products (FDP) in human blood plasma by immunoturbidimetry on coagulometers "Technology Solution". FDP s are a marker of activation of blood clotting and fibrinolysis, since they are formed in the bloodstream as a result of plasmin lysis of fibrinogen and fibrin. A high level of FDP will be registered in patients with thrombosis, thromboembolism, DIC syndrome of various genesis
A set of calibration plasmas for quantitative determination of PDF on coagulometers "Technology Solution" (TS- FDP calibrator) on the TU 21.20.23-118-42349142-2021	RZN 2022/19286 of December 30, 2022	The medical device "TS- FDP control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of FDP on coagulometers

		"Technology Solution" (TS-FDP)". The medical device "TS-FDP control" is intended only for the set "TS-FDP".
A set of control plasmas for quantitative determination of PDF on coagulometers "Technology Solution" (TS- FDP control) according to TU 21.20.23-116-42349142-2021	RZN 2022/19286 of December 30, 2022	The medical device "TS- FDP control" is intended for in-laboratory quality control of a medical device "A set of reagents for quantitative determination of FDP on coagulometers "Technology Solution" (TS-FDP)". The medical device "TS-FDP control" is intended only for the set "TS- FDP".
Disposable cuvette for the line of coagulometers "Technology Solution", in versions: - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 μl) – 1000 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K2 (700 μl) – 1500 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K3 (1 ml) - 1000 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution 190", "Technology Solution 400": K4 (1.5 ml) - 500 pcs. in the package. - Disposable cuvette for coagulometers "Technology Solution 60", "Technology Solution	RZN 2020/10843 of June 17, 2020	The medical device "Disposable cuvette for the line of coagulometers "Technology Solution"" is intended for placing blood plasma samples and reagents, mixing them and examining blood plasma samples using coagulometric analysis methods on the line of coagulometers "Technology Solution": "Technology Solution 4", "Technology Solution 60", "Technology Solution 190", "Technology Solution 400".

190", "Technology Solution 400": K 5 (2 ml) - 400 pcs. in the package.		
"Washing solution according to TU 21.20.23-001-34065548-2019 "	RZN 2019/9048 of May 04 , 2021	The medical device "Washing solution" is used as a detergent and disinfectant reagent for cleaning the dosing systems of the line of coagulometers "Technology Solution" ("Technology Solution 60", "Technology Solution 190", "Technology Solution 400"), coagulometers "AK" and "Insula 500".

When new sets of the Technology Solution line appear, users will be informed by an authorized representative in the Russian Federation of Laboratory Solutions LLC through an Official Notification, or on the official website www.labsolution.ru , www.лаб-решение .RF.

1.8.PC configuration

Control takes place under the control of PC software The recommended PC configuration is as follows:

Processor: 1.6 GHz or higher; Memory: 512 MB or higher; HDD: 160 GB or higher. Display Resolution: 1366 * 768; OS: Microsoft Windows XP or higher. The presence of at least one serial interface (used to connect the analyzer), one network interface (used to connect LIS / MIS) and two USB interfaces.

Chapter 2 Principle of Operation

2.1. Coagulation process

Figure 2-1. Scheme of the light flux of the coagulation process.

The optical system of the coagulation process uses the method of 90% optical scattering + % of the end point. The light source is an ultra bright red LED with a wavelength of 660 nm, with a high-resolution filter. The light filter directs the light flux to the sample. A photodiode receives a diffused luminous flux at an angle of 90° and converts it into an electrical signal. After converting the light flux into a digital signal, the computer processes the data, draws a coagulation curve, and determines the coagulation time using the percent end point method.

Scattering intensity


Figure 2-2. Endpoint percentage method.

Percentage is the clotting time.

At the moment when the reagent is added, and coagulation has not yet arrived, the level of scattered light will be 0%, and when the coagulation reaction is completed, it will be 100%. Based on the coagulation curve, the clotting time can be considered the time when an individual level reaches a certain test percentage. In the figure, the coagulation point is set to 50% and the user can adjust the percentage value of the coagulation point for each time, depending on the specific situation. The smaller the value, the shorter the result time will be.

Tests that are based on this method: PT, APTT, FIB, TT, PCco, PSc coagulation factors (II, V, VII, VIII, IX, X, XI, XII, XIII), etc.

2.2. Chromogenic method and immunoturbidimetry


Figure 2-3. Light flow diagram for the chromogenic method and immunoturbidimetry.

The optical system of the chromogenic method and immunoturbidimetry uses the multi-wavelength optical absorption method. The luminous flux uses the LED lamp + optical fiber mode, applying at least four wavelengths of 405/575/660/800 nm. The light from the LED lamp passes through the filter and passes through the optical fiber to the sample. A photodiode accepts light absorption and converts it into an electrical signal. After the analog-digital conversion, the computer processes the data, records the changes in optical density per minute, and uses the relationship between changes in absorption and concentration (or percentage activity) to check the concentration (or percentage activity) of the test substance.

A chromogenic test is the use of the principle of colorimetry to measure the absorption of light (405 nm) by a liquid in a cuvette. The amount of light that reaches the photo detector is converted into an electrical signal that is proportional to the activity of the enzyme.

Tests that are based on this method: FXIII, AT3, Plg, PCch, Hep, a2 AP.

The immunoturbidimetry method is based on measuring the absorption and scattering of monochromatic light by suspended particles formed during the antigen-antibody reaction. There is a change in the optical density of the medium detected photometrically. After plasma incubation for some time, a stabilizing reagent and an antigen-sensitive reagent are added. Then the reaction mixture is placed in a measuring channel with a wavelength (575/660/800 nm, depending on the parameters specified in the test protocol) and the change in light absorption caused by the antigen-antibody reaction is photometrically recorded.

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2.3. The measurement process in the presence of various interferons

The sample is tested in the presence of interferents by measuring the optical absorption of the sample at three wavelengths (405 nm, 660 nm, 800 nm). For various interferons, the methods use a certain set of wavelengths, as a rule, in each method, the main and additional wavelengths are used for the operation of the HIL function. The list of possible interferents and the correspondence of the required wavelength are presented in the table below.

Interference	Wavelength
--------------	------------

Triglycerides	660 nm
Triglycerides + hemoglobin	800 nm
Triglycerides + hemoglobin + bilirubin	405 nm

The HIL function is configured by a specialist of an Authorized Representative in the Russian Federation of Laboratory Solutions LLC during the commissioning of the analyzer.

Determination of the presence of interferons in the sample that may affect the test result is performed during sample preparation. If the operator notices that the sample has interferences, he must set a flag at the time of assignment of the task for this sample, requiring testing using the HIL function.

The purpose of testing a sample using the HILL function

To use the HILL function during sample testing, follow the steps specified in paragraph 7.5 of this manual, then set the flag in the window indicating the activation of the HIL function. In the figure below, the flag setting is highlighted in a circle.

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Figure 2-4. Activating the HILL function

Chapter 3 Installation

3.1. Unpacking the device

3.1.1. Stages unpacking

Unpack the instrument and remove material for transport. Save all packing and transportation materials in the event of a possible shipment of the instrument.

1) The box with the device should stand straight, and the arrows should be pointing up.

2) Using tools, open the side walls of the box, remove the accessories and check the equipment according to the packing list. In the event of a shortage, immediately contact the distributor or the representatives of the Laboratory Solutions Limited Liability Company. Check the contents, which should correspond to the list specified in clause 1.1.2.

3) Remove the top cover. Take the device from the packaging material by using special handles at the bottom of the device on both sides with at least 4 people.

Caution: *Accessories must match the packing list. If a component is missing, damaged, or does not match the packing list, contact your distributor or representative. Please store the packaging correctly for further packaging and shipping. The ambient transport environment for the instrument: -20 ~ 55, relative humidity up to 90%, barometric pressure 50 kPa ~ 106 kPa. During transport, pay attention to the Top, Fragile, Protect from Rain and Moisture icons. The limit of the layers on the package and transport the device in accordance with the requirements..*

3.1.2. Terms of transportation

- Before transporting the instrument, make sure that it is reset, the sampling needle is in the raised position.
- Under normal conditions, trolleys and other vehicles can be used for short distances.
- During the entire transportation of the instrument, pay special attention to the safety of the front panel, the sampling probe and the connecting elements of the fluid system on the rear panel.
- During transportation it is recommended to keep the device in an upright position, do not tilt.
- Avoid vibrations if possible. Before the subsequent operation of the device, it is recommended to check and calibrate.

Caution: *Keep the packing box for long distance transport. The device should be*

installed on a flat working surface, without corners and irregularities.

3.2. Installation conditions and operation

The automatic coagulometer "Technology Solution 190" must be installed by qualified personnel. When it is necessary to change the installation location of the device, also contact qualified specialists.

It should be noted that if malfunctions have occurred, as a result of moving the system carried out by the user, this leads to the loss of all warranty obligations, even if the warranty period has not expired.

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- Away from direct sunlight;
- In a dust-free place;
- Absence of strong electromagnetic radiation;
- Do not use this device at low relative humidity, especially if synthetic materials (synthetic clothing, carpets, etc.) are present near it, may lead to erroneous results due to the influence of electrostatic discharges.
- A grounded power outlet must be located in the immediate vicinity of the appliance;
- Leave enough space for work and maintenance. Creating conditions for heat exchange of the working device, leave approximately 50 cm from the walls to the back panel, sides, top of the device;
- The room should be well ventilated;
- Avoid excessive humidity and high temperatures;
- Avoid vibration and shaking.
- The device must be installed on a flat working surface, without inclination and irregularities.

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Attention: 1) *This device operates at an ambient temperature of 10 °C 30 °C and humidity $\leq 70\%$.*

2) *After installation, avoid moving the instrument frequently. Use a stable trolley to move it. The angle of inclination should not exceed 15 °.*

3) *Installation and relocation is carried out by qualified specialists.*

3.3. Power requirements

1) a.c. 100-240 V, 2) 50 Hz / 60 Hz, 3) 500 W

**Warning:**

- 1) The AC power supply must be well grounded (ground voltage <5V).**
- 2) The power supply must be stable when combined with high-power electrical appliances.**
- 3) If smoke, a sharp smell or unusual noise emanating from the device is detected, immediately turn off the power and contact the service center of an Authorized representative in the Russian Federation.**
- 4) To disconnect the power cord, gently pull the plug, holding directly to the cord is highly discouraged.**
- 5) Do not open the back cover, side covers and panels of the switched-on appliance.**

3.4. Installation

1) Space requirements: the distance between the right panel of the device and the wall should be ≥ 100 cm, the distance between the rear panel of the device and the wall should be ≥ 100 cm. Ensure that there is enough space on the desktop and under the table for deionized water tank and containers for waste liquid (also to provide repair and maintenance).

2) Use the WLAN interface to connect the instrument to the PC and use the electrical cable to connect the instrument to the network.

3) A container with deionized water and a container for waste liquid are connected to the interfaces on the device in accordance with the marking on the rear panel.

3.5. Calibration

Calibration of the instrument should be carried out using quality control standards and materials. The calibration test is not necessary to perform daily, but always at least once for each series, in need of calibration. Calibration data will be used until new calibration data is saved. However, due to the fact that external changes may affect the results of the tests, daily calibration is recommended to obtain more accurate test results.

Chapter 4 Launch

The list of medical devices and consumables required for use in conjunction with the medical device "Automatic coagulometer "Technology Solution 190" for its intended purpose.

To perform a study of the hemostasis system (testing) in blood plasma according to the principles specified in paragraphs 2.1 and 2.2, the device requires reagent kits, calibrators, cuvettes and a washing solution. All components necessary for the operation of the device are specified in section 1.7 of this manual.

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Disposable cuvettes (see paragraphs 1.2.4 and 1.7 of this manual)

The cuvettes are delivered in packages. K2 cuvettes are loaded into the hopper for measuring cuvettes by pouring them out of the package to the red mark located on the inside of the hopper. The cuvette can be reloaded at any time, regardless of the status of the analyzer.

The automatic coagulometer "Technology Solution 190" automatically detects the presence of cuvettes. If there is not enough cuvette in the feed line (26 cuvettes or less), a warning message is displayed on the screen. If there are no cuvettes left in the feed line, an error message appears. The analyzer stops working with new samples and performs a controlled stop. The "Warning" LED indicator on the analyzer panel turns red. In order to be able to share the cuvettes with the coagulometer and the correctness of the results obtained, it is necessary to make sure that there is a special marking on the upper face of each cuvette (see paragraph 1.2.4 of this manual).

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The washing solution is used to clean the reagent probe and the sample probe during operation. The cleaning cycle is performed according to the factory settings, and additional washing cycles can also be set in the test settings. The procedures for installing and replacing the washing solution are described in Chapter 6 of this manual.

Calibrators

A calibrator is a sample with a known analyte value obtained during production and used for calibration of the test. Depending on the specific method, sets of calibrators can be used (see paragraph 1.7 of this manual), as well as calibrators included in reagent kits. The use of calibrators is performed to obtain the calculated values of the calibrated method. Installation of calibrators on board the analyzer is carried out in the sample compartment using special cuvettes of the medical device "Disposable cuvette for the line of coagulometers "Technology Solution"", in versions:

- Disposable cuvette for coagulometers "Technology Solution 60",

"Technology Solution 190", "Technology Solution 400": K3 (1000 pcs. in the package).

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4.1. Pre-launch preparation

Before starting work on the device, the user must perform the following actions.

1. Check if there is enough distilled water.
2. Empty the waste container
3. Add the required number of measuring cuvettes
4. Empty the tank for the spent cuvettes.
5. Make sure that there are no obstacles to the movement of mechanical parts.
6. Check whether the device is connected to the network via UPS

4.2. Enabling the system

To start the device, you must first press the switch at the far wall of the device, then at the front panel, then turn on the control computer and run the device software.

4.3. Indicators

In the upper right corner of the front panel of the device there are three indicators that indicate the operating status of the device. You should pay attention to changes in the indicators during startup and during operation of the device in order to know about the operating status of the device:

- POWER: Green. When the main power switch is on, the POWER indicator lights up, it means that the device is turned on normally; if the POWER indicator is blinking, this indicates that the cooling of the reactants has not reached the desired value.
- RUNMING: Yellow. If the RUNMING indicator lights up when the secondary power switch is turned on, it means that the device is turned on and in working condition.
- WARNING: Red. When the main and secondary power switches are turned on, the WARNING indicator is not lit, which means that the instrument is operating normally; If the WARNING indicator is blinking, it means that a malfunction has occurred in the operation of the device.

4.4. Login

After turning on the device and loading the system, a window will appear to enter the system. See Figure 4-1.


Figure 4-1. Login window.

The user name is Admin, the password is 888888. After entering, you will be successfully logged in.

4.5. Main menu

After successful login you will see the main menu window. See Figure 4-2


Figure 4-2. Main menu.

The main menu is the main window for all users of the "Technology Solution 190" device, in which users can click the icon to select a specific operation.

4.6. Reference

Click the Help icon in the status bar to open the Help and Information window. See Figure 4-3.


Figure 4-3. Help and information window.

Click the Help button or the F1 shortcut key to view software help information. Click the "About Program" button to open a dialog box with information about the software. See Figure 4-4.


Figure 4-4. Information about the device software.

Chapter 5 Test Settings

Press the "Menu" button and then "System Setup" → "Tests" to enter the test setup menu. See Figure 5-1.

The screenshot shows the 'Test' menu with a list of tests on the left and configuration fields on the right. The selected test is 'LA1'. The configuration fields include:

- Test Info:** Test's Abbr (LA1), Test's Full Name (LAG-C), Test Meth (For PT), Data Format, Sensitivity.
- Ref. Range:** 0.0 to 45.0.
- Coag. pt. (%):** 50.
- Max. Time (S):** 240.
- Spl. Vol. (μL):** S1: 100, S2: 0.
- Diluent's Name:** (empty).
- Diluent's Vol. (μL):** 0.
- Pre-rinse (Cycle):** CLEAN 1 * 0.
- Post Rinse (Cycle):** CLEAN 1 * 0.
- Reagent Vol. (μL):** R1: 0, R2: 0, R3: 0.
- Reagent Name:** (empty).
- Aspirating Vol. (μL):** 0.
- Pre-rinse (Cycle):** CLEAN 1 * 0.
- Post Rinse (Cycle):** CLEAN 1 * 0.
- Other Setup:**
 - Incubation (S): 0.
 - Lower Limit: (empty).
 - Upper Limit: (empty).
 - Missing: (empty).
 - Level: M.
 - Correction Formula: E 1.000, B 0.000.

Buttons at the bottom: Save, Calculation Param, Add, Delete, Print, Restore to Default, Display Sequence.

Figure 5-1. Test setup menu.

The functions of this menu include such parameters as saving, calculating, adding, deleting, printing, restoring default parameters and choosing the sequence of information display.

Test names are displayed in the list on the left. Information about the parameters corresponding to the selected test is displayed on the right. Default test names are not subject to change..

5.1. Test parameters

Name: The name of the test.

Full name: Full name of the test.

Range of reference values: reference range of measurement results.

Coagulation point (%): The end point in the coagulation process. The default value is 50%.

Max. time: Maximum measurement time.

Start time: the start time of the measurement for the chromogenic method and immunoturbidimetry.

End time: end measurement time for the chromogenic method and immunoturbidimetry.

Sample volume S1 (μL): sample intake volume S1.

Sample volume S2 (μL): the amount of sample diluted in step 2.

Pre-rinsing: additional rinsing cycle set before sampling. If set to 0, there will be no flushing.

Post-flushing: an additional flushing cycle, established after sampling. If set to 0, there will be no flushing.

The incubation time of reagent 1: Select the interval with a minimum value of 20 seconds between the time when the sample starts to warm up and the time it takes to add reagent 1, which can be changed by an increase of 20 seconds. However, if the flush mode is selected immediately after sampling, each flush cycle will be increased by 20 seconds. In addition, if the flush mode is selected before adding reagent 1, each flush cycle will be increased by another 20 seconds.

Pre-rinsing (R1): additional rinsing cycle installed before collecting reagent 1. If set to 0, no rinsing will be performed.

Post-flushing (R1): additional flushing cycle, installed after collecting reagent 1. If set to 0, no flushing.

The incubation time of reagent 2: The minimum interval between the addition of reagent 1 and the addition of reagent 2 is 20 seconds. The interval can be extended in increments of 20 seconds. However, if flushing is set immediately after the sample is taken, for each flush cycle, the time must be increased by 20 seconds. Similarly, if the rinsing is set before adding reagent 1, each washing cycle must be increased by 20 seconds.

Pre-rinsing (R2): additional rinsing cycle, installed before collecting reagent 2. Up to 9 cycles can be set. If set to 0, there will be no flushing.

Post-flushing (R2): additional flushing cycle, installed after collecting reagent 2. You can set up to 9 cycles. If set to 0, there will be no flushing.

The incubation time of reagent 3: The minimum interval between the addition of reagent 2 and the addition of reagent 3 is 20 seconds. The interval can be extended in increments of 20 seconds. However, if flushing is set immediately after the sample is taken, for each flush cycle, the time must be increased by 20 seconds. Similarly, if the rinsing is set before adding reagent 1, each washing cycle must be increased by 20 seconds.

Pre-rinse (R3): additional rinse cycle, installed before collecting reagent 3. Up to 9 cycles can be set. If set to 0, there will be no flushing.

Post-flushing (R3): additional flushing cycle, installed after collecting reagent 3. You can set up to 9 cycles. If set to 0, there will be no flushing.

Setup analysis with re-dilution: the sample is re-tested after a given multiple dilution. Select the upper and lower limits of the result and the dilution factor. If the test results go beyond these values, the sample should be tested again after the specified multiple dilution. After selection, the system automatically starts this function.

For coagulation techniques: when the concentration is high, take a measurement after 2-fold dilution, and then multiply the result by 2; when the concentration is low, measure after a single dilution, and then multiply the result by 1/2.

For chromogenic methods and immunoturbidimetry: when the concentration is high, perform the test after 8-fold dilution, and then multiply the result by 8; when the

concentration is low, measure after 1/8 dilution, and then multiply the result by 1/8.

Determination of high / low concentration: concentration exceeding the established range is high; concentration below the set range is low.

Setting re-analysis: If you need to re-test the sample, select the upper and lower limits of the result. If the test value is outside these limits, the sample will be re-tested. After selection, the system automatically starts this function.

Change formula: Change measurement results. Formula: $Y = KX + B$; X— measurement result, Y— changeable result. If changes are not necessary, you can choose the values $K = 1$ and $B = 0$.

5.2. Saving settings and test parameters

After setting up the test, click the Save button. If the parameters are set correctly, the settings will be saved and the current menu will be updated; if the parameters are not set correctly, the settings will not be saved, a prompt window will appear with a marked error.

5.2.1. Calculation of parameters

Click the Calculate Parameters button to open the Calculate Parameters dialog box. See Figure 5-2.

Calculate	Param. Select	Param.'s Full Name	Unit	Data Format	Ref. Range	
					Lower Limit	Upper Limit
☒	PT-%	Prothrombin Activity	%	- ** -		
☒	PT-INR	PT-INR	INR	- ** -	1.0	6.0
☒	PT-Ratio	PTR	Ratio	- ** -		
☒	dFbg	Fibrinogen PT Test	mg/dL	- ** -		

Normal ISI

Figure 5-2. Parameter calculation window.

You can select the calculation parameters, change the name of the test and the full name of the test by selecting the appropriate test, you can change the data format, the lower and upper limits of the measurement range, the normal value, and the MICH value. Next, click "OK." If the parameters are set correctly, the system will save the settings, close the "Calculate Parameters" dialog box and return to the current main menu. If the parameters are not set correctly, a window will appear in which the error will be indicated. Click the "OK" button, then the error window will be closed, and the settings will not be saved; Click the Cancel button, the system will close the Calculate Parameters dialog box

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The "Calculation Parameter" window includes the possibility of entering reference ranges for each activated parameter. To interpret the test results obtained, it is necessary to refer to the instructions of the set of reagents used on the analyzer (see paragraph 1.7 of this manual).and return to the current main menu.

5.2.2. Add test

Click the Add button to open the Add dialog box. See Figure 5-3.


Figure 5-3. Dialog add test.

You can enter a test name and select a test method. Click OK. If the parameters are set correctly, the system will save the settings, close the "Add" dialog box, return to the current main menu and update the list of tests. If any parameter is set incorrectly, a

window will appear with a hint, and an error will be highlighted. Click OK, and the prompt window will close. Click the Cancel button, the system will close the Add dialog box and return to the current main menu.

5.2.3. Remove test

Click the Delete button to open the Delete dialog box. See Figure 5-4.


Figure 5-4. Dialog box delete test.

Select the test you want to remove. Press the "OK" button, the system will delete the selected test, close the "Delete" window, return to the current main menu and update the list of tests. Press the "Cancel" button, the system will close the "Delete" window without deleting the selected test, and return to the current main menu.

5.2.4. Print

Click the "Print" button to print information about the parameters of the selected test.

5.2.5. Restore Defaults

Click the Restore Defaults button to open the Restore Defaults dialog box. See Figure 5-5.


Figure 5-5. Restore default settings window.

Click "Yes", the system will restore the default settings, close the window, return to the current main menu and update the list of items. Click "No", the system will not restore the default settings, close the window and return to the current main menu.

5.2.6. Test Display List

Click the Display Test List button to open the Display Test List dialog box. See pic.

5-6.


Figure 5-6. Menu display list of tests.

Test names are displayed on the left (in the same sequence as in the test menu). Action options are displayed on the right, up, previous, next and down. Select a test and click the action button on the right. When the first item is selected, only "Next" and "Down" actions can be performed; when the last item is selected, only the "Previous" and "Up" actions can be performed.

(1) OK

Press the "OK" button, the system will change the test sequence, close the "Test Display List" dialog box, return to the current main menu and update the test list.

(2) Cancel

Click the Cancel button, the system will not change the test sequence, close the Test Display List dialog box and return to the current main menu.

Chapter 6 Reagent Settings

To carry out testing using the medical device "Automatic Coagulometer "Technology Solution 190" for its intended purpose, it is necessary to use the medical devices specified in paragraph 1.7 of this manual. The Reagent tab is used to place reagents on board the analyzer and identify them in the software. The Reagent tab displays a map of the location of reagents on the disk. Both reagents and the washing solution are placed in the reagent disk (paragraph 1.7 of this manual), each of which has a different color indication. Next to the reagent disk is a disk for setting diluents and buffers. To visually display the reagent level, a color indication of the external contour of the position is provided: blue color – normal level, purple color – low level, red color – absence of reagent.

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Click the "Reagent" button on the toolbar or the "Menu" button, and then "System Setup" → "Reagent" to enter the reagent setup menu. See Figure 6-1.


Figure 6-1. Reagent settings menu.

To install the required reagent, you need to select a position on the reagent disk, then select the name of the required reagent from the drop-down menu in the table on the right. You can enter information about reagents, including reagent rack, reagent position, reagent name, vial type, batch number, expiration date, and reagent bar code. Reagent barcode can be entered manually or scanned. The positions of reagents on the reagent tank and sample tank are shown on the left. There are 4 reagent racks in the reagent tank, and each reagent rack includes 6 reagent positions, a total of 24 reagent positions. When the position of reagent 5 is selected, the name of the reagent can only be CLEAN I or CLEAN II, and the type of vial cannot be edited. In the sample container there are 5

reagent positions, in which 1 position is for cleaning solution, you can put the reagents in other positions.

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Depending on the situation, the reagent can be replaced both during waiting and during sample testing. When the analyzer is in standby mode, it is necessary to get a bottle with a reagent that requires replacement and put a bottle with a new reagent in its place. To detect the reagent level, the analyzer will perform a check before performing the analysis automatically. If, during testing, the color indication of the contour of the reagent position has changed to purple, it means that the reagent bottle will soon be emptied. To replace the reagent, you need to press the pause button, then wait for the message that the work is suspended and replace the reagent. Then click the "Start test" button and the analyzer will continue working from the moment when it was suspended.

To expand the capabilities of the analyzer in offline mode, there is a function of automatic transition of the reagent probe from one reagent to another of the same reagent. In order for this function to be used during testing, it is necessary to install 2 or more vials of the same reagent before starting testing.

6.1. Reagent Parameters

Reagent Position: Sets the position at which the reagent will be installed.

Reagent name: reagent name.

Vial Type: Vial Type.

Series: reagent batch number.

Expiration date: Expiration date of the reagent.

Reagent barcode: reagent barcode can be entered manually or scanned.

6.2. Reagent Information

Click on the reagent designation on the reagent container or sample container, the reagent information will be displayed on the right. You can edit the information corresponding to the reagent.

6.3. Preservation

Click the Save button. If the parameters are correct, they will be saved; if any parameter is set incorrectly, a window will appear with a hint, and the first error will be highlighted, and the reagent parameters will not be set. Click "OK" to close the window with a hint.

6.4. Barcode Scanning (Optional)

To identify and safely combine medical devices (reagents), it is required to scan the barcode applied to each reagent vial. Click the "Barcode Scanning" button, the "Barcode

Scanning" dialog box corresponding to the selected reagent opens. See Figure 6-2.


Figure 6-2. Barcode scanning window.

You can choose the barcode scanning range: all reagent positions or preset positions. Click the Start button and the system will load the barcode to the specified reagent, close the Barcode Scanning dialog box and return to the current main menu. Click the Cancel button, the system will not scan the barcode, close the Scan dialog box and return to the current main menu.

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The barcode system of the device is compatible only with the line of kits and reagents "Technology Solution".

6.5. Deletion

Click the Delete button to open the Delete dialog box. See Figure 6-3.


Figure 6-3. Delete dialog box.

You can use the Shift and Ctrl keys on your keyboard to select multiple reagents. Click OK and the system will delete the selected entry, close the Delete dialog box and return to the current main menu. Click the Cancel button, the system will not delete the selected entry, close the Delete dialog box and return to the current main menu.

Chapter 7 Testing procedure

Before starting the testing, it is necessary to perform the preparatory work described in paragraph 4.1 of this manual.

Click the Test button on the toolbar or the Menu button, and then the Test option to enter the test menu. See Figure 7-1.


Figure 7-1. Testing menu.

The main features provided in this menu include sample selection, counting, printing, transfer, bar code scanning (optional), retest, response curve, etc. The functions provided in this menu vary depending on the level of user access.

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1. Prepare the primary test tubes by removing the closing caps to avoid collision and possible damage to the sample probe.
2. Place the prepared test tubes in the sample compartment. If the test tubes have a barcode, then it is necessary to expand the barcode towards the slot in the sample disk.
3. In the software interface, specify the necessary testing methods for each sample. If the analyzer is connected to the LIS, then it is necessary to scan barcodes. After successful barcode scanning, the analyzer will receive the corresponding task from the LIS database.

4. After all samples have been installed and methods have been obtained or set for them, click the "Start test" button located in the upper right corner of the tab.

During the testing process, the analyzer will inform the operator about the status of the samples being performed, as well as display the completed test results. After the analyzer finishes the testing procedure, a message will appear in the lower status bar. The procedure for calculating the displayed results is described in Chapter 2 of this manual. To interpret the results obtained, refer to the instructions of the reagent kits (see paragraph 1.7 of this manual).

For the rules for performing repeated measurements, see paragraph 7.5.9 of this manual.

7.1. Sample container

Here the sample number is displayed with different colors indicating different sample types. The color of the sample indicates its test state. The selected sample position in the sample container, the sample position in the sample list will be selected accordingly, and information on the corresponding test will be displayed in the list of details.

7.2. STAT-test

A solid color indicates whether there is a sample or not. When the sample position in the tank for STAT samples is selected, the corresponding information in the list of samples will be selected, and information about the test corresponding to the sample will be displayed in the list of details.

7.3. Sample list

Here you can enter information for a sample or QC (quality control), information on test results and QC measurements.

7.4. Test Information

Information about the parameters of the selected test is displayed here. If several tests are selected, the parameters of the last selected test are displayed.

7.5. Sample input

Click the Sample Input button to open the Sample Input dialog box. See pic. 7-2.

Figure 7-2. Sample entry menu.

In this menu, you can specify the following information: sample position, sample number, sample type, test tube type, name, medical history number, gender, etc. In the "Tests" section, you can select tests to be performed for the sample.

Sample Position: The position corresponding to the sample in the sample container, positions 1-55.

Sample Number: Used to identify the sample.

Sample type: sample type.

Tube Type: The type of tube containing the sample, including the conventional tube or cuvette.

Name: the name of the patient.

Case Number: The patient's case number.

Gender: Gender of the patient.

Date of birth: date of birth of the patient.

Age: age of the patient.

Patient type: patient type, System settings -> Data dictionary -> Patient type setup.

Section: System settings -> Data dictionary -> Section settings.

Hospital name: System settings -> Data dictionary -> Hospital settings.

Bed number: The bed number can only be set up if the patient type is "hospitalized".

Time of submission: the time when the sample is sent to the study.

Testing time: sample testing time.

Customer: the operator who provides the sample for testing.

Artist: An operator who is testing a sample.

Verifier: the operator that checks the sample.

Diagnosis: information about the diagnosis of the sample.

Remarks: Note to the sample.

Test Information: All tests provided by the system are displayed on the left. The test is displayed as a button. Tests that are ready to run are highlighted; tests that are not ready

are displayed in gray, and the reason why they cannot be performed will be displayed when you hover over them. The name of the test profile is displayed on the right in the form of a button. Click the button and the tests contained in the profile will be automatically selected; pressing it again will cancel. If the tests contained in the test profile cannot be executed, the tests will not be selected..

7.5.1 Preservation

Click the Save button. If the sample data is set correctly, the settings will be saved. The number of default analyzes is one. If any sample information is incorrectly installed, a window will appear with a hint, and the first error will be highlighted.

7.5.2. Batch input

Click the Batch Input button to open the Batch Input dialog box. See Figure 7-3.


Figure 7-3. Batch input menu.

Select the start and end position of the samples, click "OK". If the starting position of the sample is greater than the ending position, a pop-up message will appear; If no test sample is selected, a pop-up message will appear. When all settings are correct, enter sample information between the initial position of the sample and the final position..

7.5.3. The number of analyzes

Click the "Number of analyzes" button to open the dialog box. See Figure 7-4.


Figure 7-4. The "Number of analyzes" dialog box.

Click in the field "Number of analyzes" to edit the number of analyzes of the corresponding test (1-10). Select the field "Number of analyzes" for all tests, then the number of analyzes of all tests will be set the same.

Press the "OK" button, the system will close the "Number of analyzes" dialog box, saving the number of analyzes, return to the current main menu and update the list of samples and the list of tests. Click "Cancel", the system will close the "Number of analyzes" dialog box, without changing the number of test analyzes.

7.5.4. Report

To check the sample that passed the measurement you need to click the "Report" button. A successfully measured sample will be marked with a "✓" sign and there will be no possibility of changing it and re-analyzing it.

7.5.5. Print

Click the Print button to open the Print dialog. See Figure 7-5.


Figure 7-5. Menu print measurement results.

After setting the print range and format, click OK. If the print settings are correct, the system will close the "Print" dialog box, print the record and return to the current main menu. If you click Cancel, the system closes the Print dialog box and does not print the entry, returning to the current main menu.

7.5.6. Data transfer to LIS

If the system is not connected to the LIS server, then when you click the Transfer button, a dialog box will appear asking you. See Figure 7-6.


Figure 7-6. Connection to the LIS server.

After successfully connecting to LIS, the "Transfer" dialog will appear. See pic. 7-7.

Figure 7-7 Sample Data Transfer Window

This window sets the range of records to be transferred. If the connection is established, click the Start button. If you selected a record in which there is no measurement data for this sample, a dialog box appears with the message "There is no data record in the selected range!". After clicking the "OK" button, the dialog box will be closed and the "Transfer" window will open. If the selection and input is correct, the system will begin to transfer the selected entries, close the Transfer dialog box and return to the current main menu. Pressing the "Cancel" button, the system will not perform data exchange, close the "Transfer" dialog box and return to the current main menu.

7.5.7. Barcode scanning (optional)

Click the "Barcode Scanning" button to open the "Barcode Scanning" dialog box. See pic. 7-8.


Figure 7-8. Barcode scanning window samples.

This window sets the range of bar codes to be scanned. After selecting a range, click the Start button. If all settings are correct, barcode scanning will begin. During the scan, the status bar displays "Scan Barcode ...". After the scan is complete, the message will not appear in the status area. After scanning the bar codes, the scanned information will be entered into the appropriate edit field and automatically saved. For a sample that does not have a barcode, a default value will be displayed in the edit box. Clicking Cancel will close the Barcode Scan dialog box.

7.5.8. Skip positions

Click the "Skip Position" button to open the "Skip Positions" dialog box. See Figure 7-9.


Figure 7-9. Window "skip positions."

In this window you can select "Skip the selected sample position" or "Skip all positions". If you selected "Skip the selected sample position", the sample in the selected position will be skipped (its status will be: not tested or not completed); If you selected "Skip all positions", the system will skip samples in all positions (sample status will be: not tested or not completed), close the Skip Positions dialog box and return to the main menu. When you click Cancel, the system will close the "Skip Positions" dialog box and return to the main menu.

7.5.9. Retest

Samples that are not fully tested cannot be retested. When the device is in pause mode, the repeated test is not allowed.

Click the "Repeat Test" button to open the "Repeat Test" dialog box. See Figure 7-10.


Figure 7-10 Retesting window.

In this window, you can select items for retesting, as well as ways to process the results of retesting. In case of obtaining a result that raises doubts about the correctness, including if the influence of interferences is suspected, a confirmation test should be performed. To do this, you need to activate the HILL function when performing a repeated measurement. To perform repeated testing, click OK. When you click on the "Cancel" button, the system will close the "Repeat" dialog box and return to the main menu.

7.5.10. Reaction curve

Click the Reaction Curve button to view a graph of the reaction curve corresponding to the selected test.

Chapter 8 Quality Control

8.1. Quality Control Settings

The automatic coagulometer "Technology Solution 190" can display 365 quality control points. With the possibility of simultaneous storage of 24 types of control materials.

Press the Menu button → Quality Control → Quality Control Settings to enter the quality control settings. See Figure 8-1.


Figure 8-1. Quality control setup menu.

The names of the tests are displayed on the left. Information about Quality Control and information about the parameters corresponding to each test are displayed on the right. If the quality control parameters are set for the selected test, they will be displayed in the parameters; if the parameters are not set, the default values are displayed. The entire list consists of 24 variants of Quality Control materials.

The evaluation of the control results is carried out using the Westgard control rules. For this purpose, one measurement is carried out daily for each of the 2 control materials (with normal and pathological content of indicators) or two measurements in the same control material. In the process of evaluating the results, Levy-Jennings graphs are used for each of the control materials.

8.2. Quality control materials

Click the Quality Control Materials button to open the Quality Control Materials dialog box. See Figure 8-2.

Figure 8-2. The window of quality control materials.

QC is the selected file from the quality control list. You can set the QC name of the control material, the series, the level and shelf life. Click OK. If all the settings are correct, the system will close the QC dialog box, save the settings of the QC parameters and return to the main menu. If any parameter is incorrect, a window will appear with a hint, and the first error will be highlighted. Press the "Cancel" button, the system will close the quality control dialog box and return to the main menu.

8.3. Preservation

When all the parameters and information are set correctly, click "OK". If all settings are correct, the system will save the settings and update the current quality control list. If any parameter is incorrect, the corresponding prompt window will appear with the first error marked, and the settings will not be saved.

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8.5. Quality Control Charts (Levi Jennings Charts)

Click the "Graph. KK" on the toolbar or the "Menu" button → KK → KK Chart to enter the menu of quality control charts. See Figure 8-3.


Figure 8-3. Quality control chart menu

On the left are options for quality controls and information about quality control parameters, including the name of the test, the QC file, the level and expiration date. Quality control graphs (Levy-Jennings graphs) are displayed on the right. The KK table shows the number of KK points (N) and the calculation of averages, SD (mean square deviation) and CV (%) (coefficient of variation) of Levi-Jennings maps, on which the measurement results of control materials are applied daily or from series to series. Such graphs are constructed for each test and for each level of control material.

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Figure 9-1. Report menu.

In the interface, you can enter / edit patient information, view and check test results, print out test reports and display historical data of all samples and STAT samples.

Chapter 10 Calibration

10.1 Calibration curve

Press the "Menu" → "Calibration" → "Calibration Curve" button to enter the "Calibration Curve" menu. See Figure 10-1.


Figure 10-1. Menu calibration curve.

You can view calibration information that matches the selected test, including some calibration settings, calibration data, and a calibration curve. Click the "Print" button to print the calibration data and the calibration curve of the parameter displayed on the current interface..

10.2 Automatic calibration

Press the "Menu" → "Calibration" → "Automatic Calibration" button to enter the automatic calibration menu. See Figure 10-2.


Figure 10-2. Auto Calibration Menu.

10.3 Settings

Click the Setup button to enter the dilution factor and default value. See Figure 10-3.

The 'Setup-D-Dimer' window contains the following settings:

- Calibrator: Dimer
- Batch No.: 234234
- Validity: 7/31/2015
- Dilution Mode: Auto
- Dilution Group: 1:1/1, 1/2, 1/4, 1/8, 1/16
- Std. Substance Value: 100.0
- Reagent WS: WS 1
- Diluents: IBS

S/N	Dilution Coefficient	Std. Substance Value	No. of Times
1	1/1	100.0	1
2	1/2	50.0	1
3	1/4	25.0	1
4	1/8	12.5	1
5	1/16	6.3	1
6			1

Figure 10-3. Menu for setting the automatic calibration mode.

You can select the dilution mode and the group of dilution factors, enter the value of the standard of the substance. When the dilution mode is set to "Manual", the second field with the list of dilution factors will not appear. When the dilution mode is set to "Automatic," the parameters of the dilution factor group can be selected from the second

drop-down list. The parameters of the group of dilution factors are divided into the parameter (user-defined) and the preset dilution parameters from the drop-down list. Select the diluent and calibrator from the appropriate lists, enter the calibrator batch number and its standard value. If you need to enter a group of calibrators with already known parameters (not requiring automatic dilution), the dilution mode is set to "Manual". Then the values of each calibrator are entered.

Click OK. If the data entry is correct, the system will save the settings, close the installation dialog box and return to the current main interface. If the input is not valid, a prompt will appear.

When the calibration setting is complete, click the "Start Test" button, enter the sample location information in the calibration standard in the "Start" dialog box and start the calibration. The test status and results will be displayed in the data list on the calibration interface. If there are several tests, the results of each test will be displayed in the list.

Click the "Reaction Curve" button to view the reaction curve of each control point. When the test is completed, the calibration curve will be automatically drawn in accordance with the selected calculation method.

Click the Save Curve button to save the current calibration curve and update the calibration curve view as confirmed.

Click the "Print" button to print the calibration data and curve.

10.4 Manual calibration

Press the "Menu" → "Calibration" → "Manual Calibration" button to enter the manual calibration settings menu. See Figure 10-4.


Figure 10-4. Manual calibration menu

After selecting an item, you can select parameters, enter calibration data, and select calculation methods. Click the Curve button. If the calibration data entry is valid, the calibration curve will be displayed; otherwise, a pop-up message will appear. Click the

Save Curve button to save the calibration curve and update the Calibration Curve interface as a valid calibration curve. Click the "Print" button to print the calibration data and curve.

10.5 Calibration history

Press the Menu button → Calibration → Calibration History to enter the Calibration History menu. See Figure 10-5.


Figure 10-5. Menu history calibration.

Different calibration histories can be viewed depending on user privileges. The administrator can view all calibration histories performed by the administrator after the instrument leaves the factory; A user service engineer and above can view all calibration histories.

Click the "Print" button to print all the entries in the calibration history list, but the calibration data will not be printed.

Chapter 11 Analytical characteristics

11.1 Coefficient of variation (reproducibility).

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11.4 Influence of interferences.

The measurement process in the presence of various interferences is specified in clause 2.3. Also refer to the operating instructions for the corresponding set of reagents, including those specified in clause 1.7 of this manual.

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Chapter 12 System Settings

12.1 Additional features

Press the button "Menu" → "System Setup" → "Additional Functions" to enter the configuration menu of additional functions. See pic. 11-1.


Figure 12-1. Setup menu for additional functions.

You can set the sample number input mode (including manual input and automatic increase in quantity). For auto zoom mode, you can also choose whether to use the prefix for numbering or not; "After checking one record, automatically go to the next record to be audited"; Re-check; Alert; Sample / Patient Information; time before hibernation; does the halogen lamp turn off during sleep mode, etc.

Click "Save", and the system will determine the correctness of the settings. If all settings are correct, the settings will be saved; if any parameter is incorrect, a pop-up message will appear with the first error noted.

12.2 Print

Press the "Menu" → "System Setup" → "Print" button to enter the print settings menu. See Figure 11-2.


Figure 12-2. Print setup menu.

You can choose the default printer, report title, report format, reports, and automatic printing (automatic on / off / printing after auditing), etc.

Click the Save button. If all settings are correct, print settings will be saved; if any parameter is incorrect, a pop-up message appears with the first error noted, and the print settings will not be saved.

Click the "User Defined" button to open the custom print content setting window. See Figure 11-3.


Figure 12-3. User-defined interface for printing content.

You can set the print content, which fields will be printed in the header and which

will be printed on the footer, and also adjust the printing order.

12.3 User Defined Parameters

Click "Menu" → "System Setup" → "User-Defined Parameters" to open the settings window for entering user parameters. See Figure 12-4.


Figure 12-4. A menu of user-defined options.

When there is no data in the list, only the "Add" button will be highlighted, and the "Change" and "Delete" buttons will be highlighted in gray.

Click the Add button to open the Add dialog box. See Figure 12-5.

Figure 12-5. Menu to add custom settings.

You can set the full parameter name, abbreviated parameter name, result variants, units, data format, reference value, and range of reference values.

Click OK. If all settings are correct, the system will close the Add dialog box, save the user-defined parameter values, return to the current main interface, and update the list of user-defined parameters. If any parameter is incorrect, a window will appear with the request, and the first error will be highlighted. Click Cancel, and the system closes the Add dialog box, does not add a custom parameter, and returns to the current main interface.

Click the Edit button to open the Edit dialog. See Figure 12-6.

Edit

Param.'s Full Name: ALB

Param.'s Abbr.: ALB123

Options of Result: 1,2,3,4,5,6,7,8
(If there is more than one option, separate with ",,")

Unit: mg/dl

Data Format: ***

Ref. Value: 1.231

Ref. Range: 1.000 - 6.000

OK Cancel

Figure 12-6. Editing user interface settings.

Information about the various parameters displayed in the Edit dialog box is the parameter of the selected user parameter. If the parameter has not been set, the default value that can be changed by the user will be displayed.

Click OK. If the settings for all parameters are correct, the system closes the Edit dialog box, saves the parameter change and returns to the current main interface. If any parameter setting is incorrect, a prompt window will appear with the first error noted.

Click Cancel, and the system closes the Edit dialog box, does not change the settings, and returns to the current main interface.

Click the Delete button, and the corresponding dialog box will appear. Click the "OK" button, and the system deletes the parameter selected by the user, closes the prompt window, returns to the current main interface, and updates the list of user-defined

parameters. Click the Cancel button and the system will not delete and close the invitation window and return to the current main interface..

12.4 Profile installation

Click "Menu" → "System Settings" → "Profile Settings" to enter the profile settings interface. See Figure 12-7.


Figure 12-7. Interface installation profile.

You can add, edit and delete profiles.

When there is no profile in the list, the "Add" button is highlighted, and the "Change" and "Delete" buttons are not available.

12.4.1 To add

Click the Add button to open the Add dialog box. See pic. 12-8.


Figure 12-8. Menu add profile.

You can set the name of the profile and select the elements that should be contained in the profile. Click OK. If all settings are correct, the system closes the Add Profile dialog box, saves the settings, returns to the current main interface, and updates the list of profiles. If any parameter is incorrect, a window will appear with the request, and the first error will be highlighted. Click Cancel, and the system closes the Add dialog box, does not add a profile, and returns to the current main interface..

12.4.2 Editing

Click the Edit button to open the Edit dialog. See Figure 12-9.


Figure 12-9. Profile editing menu.

The information displayed in the Edit dialog box is the parameters set for the selected profile. You can change the profile name and the elements contained in the profile. Click OK. If all settings are correct, the system will close the "Edit" dialog box, change the selected profile, return to the current main interface and update the list of profiles. If any parameter is incorrect, a window will appear with the request, and the first error will be highlighted. Click the Cancel button and the system will close the Edit dialog box, will not change the profile and return to the current main interface.

12.4.3 Deletion

Click the Delete button to open the Delete Profile dialog box. See Figure 12-10.


Figure 12-10. Profile removal interface.

Click "Yes" and the system will delete the selected profile, close the invitation window, return to the current main interface and update the list of profiles. Click the "No" button and the system will not delete the selected profile and close the invitation window and return to the current main interface.

12.5 Data Dictionary

Press the "Menu" button → "Settings" → "Data Dictionary" to enter the data dictionary interface. See Figure 12-11.


Figure 12-11. Data dictionary menu.

12.5.1 Adding data to the dictionary

Select the appropriate data dictionary and click the Add button to open the Add dialog box. See Figure 12-12.

Figure 12-12. Interface for adding a data dictionary.

Enter a name, a shortcut code and notes and check the box "When entering a sample, this element will be set by default". Click OK. If all the settings are correct, the system will perform the addition, close the Add dialog box, return to the current

main interface and update the list. If any parameter is incorrect, a window with a request will appear, and the first error will be highlighted. Click the Cancel button, and the system will close the Add dialog box, will not perform the addition, and will return to the current main interface.

12.5.2 Deleting data from the dictionary

The default value cannot be deleted, and only the added content can be deleted. Click the "Delete" button to open the "Hint" window. See Figure 11-14.


Figure 12-14. The menu for deleting an entry in the data dictionary.

Click the "Yes" button and the system will perform the deletion, close the "Request" dialog box, return to the current main interface and update the list. Click the "No" button, and the system will not perform the deletion and will close the "Request" dialog box and return to the current main interface.

12.6 Connection

Press the "Menu" → "System Setup" → "Communication" button to enter the communication menu. See Figure 12-15.


Figure 12-15. Communication menu.

Configuring the interface includes configuring the IP address and setting the port. The parameters in the communication mode include "Communication method", "Send test results in real time", "Wait for the response to the communication", "Communication timeout" and "Connect a workstation at startup".

12.6.1 Saving communication settings

Click the Save button. If all settings are correct, connection settings will be saved. If any parameter is incorrect, a window will appear with a hint, the first error will be highlighted and the settings will not be saved.

12.6.2 Connection

Click the "Connect" button to automatically connect to the LIS server. When the connection is successful, the LIS connection icon in the status bar will be updated; when the connection fails, the LIS connection failure will be displayed in the instrument / status test area.

12.6.3 Disconnection

Click the "Disconnect" button to automatically disconnect the connection to the LIS server. When the disconnect is successful, the LIS connection icon in the status bar will be updated, and the LIS connection failure will be displayed in the tool / status test area.

12.7 Barcode Information

Press the Menu button → System Setup → Barcode Information to enter the barcode information interface. See pic. 12-16.


Figure 12-16. Barcode information interface.

You can set the sample barcode and reagent barcode. Barcode types are displayed on the left (sample barcode and reagent barcode). The right side shows the parameter settings corresponding to the barcode. After setting the parameters, press the "Save" button. If all settings are correct, the settings will be saved. If any parameter is incorrect, a window will appear with a hint, the first error will be highlighted and the settings will not be saved.

12.8 User management

Click "Menu" → "System Setup" → "User Management" to enter the user management interface. See Figure 12-17.


Figure 12-17. Manage users menu.

The user management interface provides various functions depending on user rights. Ordinary users do not have the "Add", "Edit" or "Delete" function, and only a specific user who has registered in the system has the "Change Password" function. Administrator and above have the functions of adding, editing, deleting and changing the password.

12.8.1 Add user

Click the Add button to open the Add dialog box. See Figure 12-18.


The image shows a Windows-style dialog box titled "Add". It contains several input fields and controls. On the left, there are labels for "User Name", "Name", "Pwd.", "Confirm Pwd.", "Privilege", and "Remarks". To the right of these labels are the corresponding input areas: a single-line text box for User Name, a single-line text box for Name, a single-line text box for Pwd., a single-line text box for Confirm Pwd., a dropdown menu for Privilege (currently showing "Normal User"), and a multi-line text box for Remarks. To the right of the Privilege dropdown is a checkbox labeled "Open audit authority". At the bottom of the dialog box are two buttons: "OK" and "Cancel".

Figure 12-18. Interface add user.

You can set the username, name, password, rights, whether to open the rights of verification, comments, etc. After installation, click "OK". If all settings are correct, the system closes the Add dialog box, adds the user, returns to the current main interface, and updates the list. If any parameter is incorrect, a dialog box will appear asking you, the first error will be highlighted and the user will not be added. Click Cancel, and the system closes the Add dialog box, does not add a user, and returns to the current main interface.

12.8.2 Editing user

Click the Edit button to open the Edit dialog. See Figure 12-19.


Figure 12-19. User editing interface.

You can change the username, name, rights and notes. After making changes, click OK. If all settings are correct, the system closes the Edit dialog box, saves the edited information, returns to the current main interface, and updates the list. If any parameter is incorrect, a dialog box will be displayed with a hint, the first error will be highlighted and no changes will be made. Click the Cancel button and the system closes the Edit dialog box, does not save the edited information, and returns to the current main interface.

12.8.3 Deleting user

Click the Delete button to open the Request to Delete User dialog box. See pic. 12-20.


Figure 12-20. The delete request window.

Click "Yes", and the system will delete the user, close the "Query" dialog box, return to the current main interface and update the list of users. Click "No", and the system will not delete the user and close the "Query" dialog box and return to the current main interface.

12.8.4 Change Password

Click the Change Password button to open the Change Password dialog box. See Figure 12-21.


A screenshot of a Windows-style dialog box titled "Change Pwd.". It contains three text input fields labeled "Old Pwd.", "New Pwd.", and "Confirm New Pwd.". At the bottom, there are two buttons: "OK" and "Cancel".

Figure 12-21. Interface windows change password.

Enter the parameters and click OK. If all settings are correct, the system will close the password change dialog box, change the password, return to the current main interface and update the list of users. If any parameter is incorrect, a dialog box will appear asking you, the first error will be highlighted and the password will not be changed. Click Cancel, and the system closes the Change Password dialog box, does not change the password, and returns to the current main interface.

12.9 Lab Information

Click "Menu" → "System Setup" → "Laboratory Information" to enter the laboratory information settings interface. See Figure 12-22.

The screenshot shows a software interface for entering laboratory information. It includes a menu bar at the top with options like 'Lab Info', 'Test', 'Reagent', 'Report', 'QC Chart', 'Main', and 'Exit'. The main area contains several input fields for hospital and lab details, instrument specifications, and service contacts. A 'Save' button is located at the bottom center of the form.

Figure 12-22. Laboratory Information Interface.

You can set the name of the hospital, laboratory name, person in charge, contacts, postal code, service contact person, customer service contacts, notes, device model, instrument serial number, and installation date. After setting the information parameters, click the Save button. If all parameters are set correctly, the parameters of the laboratory information will be saved. If any parameter is set incorrectly, a window will appear asking you, the first error will be highlighted..

12.10 Testing method

Only a customer support engineer and above has the Test Method feature. Press "Menu" → System Setup → Test Met. to enter the test method interface. See pic. 12-23.

The screenshot displays the 'Test Method' configuration screen. It features a sidebar on the left with buttons for different test methods: 'For PT', 'For APTT', 'For TT', 'For Fibrinogen', 'For AT-3', 'For D-Dimer', and 'For FDP'. The main area is a table of parameters for the selected method, with columns for 'Parameter Name', 'Type', 'Value', and 'Unit'. Parameters include 'Coagulation Curve Jumping Threshold', 'Coagulation Curve Saturation Threshold', 'Coagulation Curve Two-Point Threshold', 'Differential Width', 'Differential Smoothing', 'Max. No. of Spk. (Spk/10.15sec)', and 'Spk. Rising Time (s)'. At the bottom, there are buttons for 'Save', 'Add', 'Delete', and 'Restore to Default'.

Figure 12-23. Test method interface

The main functions provided on the interface include "Add", "Delete" and "Change"; the main parameters of the test methods.

12.10.1 Saving test method settings

Select a test method. If parameters were specified for the test method, the parameter set will be displayed on the right, otherwise the default value will be displayed. You can change the parameters. After setting the parameters, press the "Save" button. If all parameters are set correctly, the settings will be saved; otherwise the prompt window will appear with the first error highlighted, and the settings will not be saved.

12.10.2 Adding a test method

Click the Add button to open the Add Test Method dialog box. See Figure 12-24.


Figure 12-24. Interface add test method.

Enter the test method, Select a classification (coagulation, chromogenic or immunological methods) and click OK. The system adds the test method, closes the Add dialog box, returns to the current main interface, and updates the interface. Click Cancel, and the verification method will not be added, and the Add dialog will close.

12.10.3 Delete Test Method

Select the verification method on the left and click the "Delete" button to open the "Request to delete a method" window. See Figure 12-25.


Figure 12-25. Request window to delete test method.

Click the "OK" button and the system will remove the test method, close the query window, return to the current main interface and update the interface. Click the Cancel button and the system will not delete the test method and close the invitation window and return to the current main interface.

Chapter 13 Service

13.1 Maintenance

Click the "Maintenance" button on the toolbar or the "Menu" → Tools → Maintenance button to enter the system maintenance interface. See Figure 13-1.


Figure 13-1. System maintenance interface

You can perform the following operations (by clicking the appropriate button):

Sleep mode, wake up from sleep mode, flush sample probe, flush reagent probe, pump fluid system, flush fluid system, reload sample probe, reagent probe reboot, reboot capture unit, reagent tank reload, reload sample tank, reload STAT tank, reload mechanical parts of the device and reboot the entire analyzer.

13.1.1 Sleeping mode

Click the Sleep button to perform the Sleep operation.

13.1.2 Exit Sleep Mode

Click the Exit Sleep button to complete the wake-up sleep operation.

13.1.3 Washing the sample probe

Click the "Flush sample probe" button to perform a flush operation on the sample probe.

13.1.4 Washing the reagent probe

Click the "Flush Reagent Probe" button to perform the flush operation of the reagent probe.

13.1.5 Filling the liquid system

Press the "Long flush" button. Check if there is enough water in the distilled water tank. If there is not enough, the system will suggest "Adding a sufficient amount of distilled water to the container before performing a long rinse." If there is enough water in the tank, a long-term flushing will be performed.

13.1.6 Emptying the liquid system

This feature is only available to users with administrator access level and higher. Click on the "Emptying the liquid system" button to clean the liquid system.

13.1.7 Reloading the sample probe

Click the "Restart sample probe" button to restart the sample probe, a progress bar appears. When the reboot is complete, the progress bar will automatically close.

13.1.8 Reloading the reagent probe

Click the reboot probe reboot button to reboot the reagent probe, a progress indicator window will appear. When the reboot is complete, the progress bar will automatically close.

13.1.9 Reload Capture Block

Click the Reboot Capture Block button to reload the capture block, and a progress bar will appear. When the sequence is completed, the progress bar will automatically close.

13.1.10 Reload reagent tray

Click the Reboot Tray Reset button to reboot the reagent tray, and a progress bar appears. When the sequence is completed, the progress bar will automatically close.

13.1.11 Sample Capacity Reload

Click the "Reboot Sample Capacity" button to reload the sample capacity, and a progress bar will appear. When the sequence is completed, the progress bar will automatically close.

13.1.12 Reloading STAT Capacity

Click the "Reboot STAT Capacity" button to reset the STAT Capacity, and a progress bar will appear. When the sequence is completed, the progress bar will

automatically close.

13.1.13 Reloading the mechanical part of the device

Click the "Restart Mechanical Part of the Instrument" button to reset the mechanical part of the instrument, and the progress bar appears. When the sequence is completed, the progress bar will automatically close.

13.1.14 Reboot the device

Click the "Restart the entire analyzer" button to restart the entire analyzer, and the progress bar will appear. When the sequence is completed, the progress bar will automatically close.

13.2 Data Maintenance

Press the Menu button → Tools → Data Maintenance, to enter the data maintenance interface. See Figure 13-2.


Figure 13-2. Data service interface

The functions provided on the interface include: backup and archive of standard samples, backup and restore of the database, backup and restore of configuration parameters, cleanup of the archive of data, QC data, log and all data, export of standard reference data, sample data QC, log, calibration information, status information and commissioning information.

13.3 Software update

The automatic coagulometer "Technology Solution 190" uses software version V1.0r. If a new version of the software is released, the update is carried out by qualified specialists of the authorized representative in the Russian Federation of Laboratory Solutions LLC.

13.4 Journal

Click "Menu" → Tools → Log to enter the log interface. See Figure 13-4.


Figure 13-4. Log interface

Log types are displayed on the left. The right is the contents of the log. Different types of logs are displayed depending on the level of user access rights. Fault information is visible only to the user at the administrator level and above. The execution of the sequence is seen only by the user of a user service level engineer and above. In the same journal type, high-level users can view lower-level user logs, but not vice versa. Logs are displayed in decreasing order of time / date. You can query all the journals or journals in the specified date range. Click the Refresh button. If there are logs, the logs will be displayed in the list; if there is no journal, the list will be empty.

13.5 Maintenance

This feature only works for a user with the rights of a customer support engineer and above. Press the Menu button → Tools c Maintenance to enter the maintenance interface. See Figure 13-5.

Figure 13-5. Maintenance Software Interface.

Chapter 14 Analyzer Service

14.1. Maintenance of the device performed by the user

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14.1.1. The procedure for daily maintenance

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- Turn off the power of the device.
- Empty the container with the used measuring cuvettes.
- Close the reagent bottles and put them in the refrigerator.
- Empty the waste liquid container.
- Fill the distilled water tank.
- Wipe the sample/reagent probes with an alcohol-free lint-free cloth.
- Wipe the working surfaces of the appliance.

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14.1.2. The procedure for conducting monthly maintenance

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14.1.3. Maintenance when the analyzer is idle

The downtime of the analyzer is considered to be the absence of any of its activities (sample analysis, calibration, etc.) for more than 7 days. After this period, it is necessary to turn on the analyzer and perform the "Long-term flushing" procedure. At the same time, it is necessary to make sure that the distilled water tank is filled, and open vials with washing liquid are installed in the reagent stand and in the washing zone of the sample probe. In case of further continuation of downtime, it is necessary to carry out the above procedure every 7 days.

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The operation "Long-term flushing" is described in paragraph (12.1.5) of this manual.

14.2. Maintenance of the device carried out by a specialist

This model of the device needs to be serviced every 6 months. Certified specialists who have been trained at the manufacturer's factory and have a valid certificate of admission to this model of the device have access to this type of maintenance.

Chapter 15 System Failures and Remedies

This chapter describes common malfunctions and how to troubleshoot them. If the problem persists according to the recommendations in this chapter, or for more information, contact Customer Service.

Malfunction	Solutions
1) The analyzer does not turn on.	Check if the plug is plugged in. Check fuse. Check the voltage.
2) The lamp in the test area is off.	Check the power supply before replacing the lamp. If the computer and instrument can be started up normally, replace the lamp.
3) Programs can not be downloaded.	Turn off the instrument and restart it after 10 seconds.
4) The printer does not turn on.	Check if the plug is plugged in. Check the "ON / OFF" button. Check fuse.
5) The printer does not print.	Check if there is a connection to the computer.
6) There is no liquid in the cleaning system.	Check if the tank is closed for cleaning correctly. Check if the high pressure pump is operating normally. Check the infusion tube.
7) The liquid in the cleaning system is full.	Check if the gravity valve block is blocked Check for leaks in the tube.
8) No data from cell detection.	Check whether the lamp is lit in the test zone. Check if there is an output signal from the preamplifier board. Check the settings of the light flux settings.
9) The magnitude of the optical signal is too high and saturation occurs.	Check the settings of the light flux settings. Check if there is an item placed in the test area. Check if the aspiration needle is blocked.

- | | |
|---|---|
| 10) The reproducibility of the results is unsatisfactory. | <p>Used contaminated reaction cuvette.</p> <p>The reaction liquid volume is less than 150 μl</p> <p>Check if the reagent is appropriate.</p> <p>Check if the sample matches.</p> <p>Replace the lamp.</p> <p>Bad contact of the signal cable.</p> <p>The reaction liquid is contaminated.</p> <p>Check if the reagent has not expired.</p> <p>Check whether the detection cell is contaminated and there is no waste liquid or residues in it.</p> |
| 11) The suction volume of the sampling probe is not constant. | <p>Check if the aspiration needle is blocked.</p> <p>You may need to replace the handset.</p> |
| 12) QC is not within specified limits. | <p>Check the expiration date of the QC fluid and that the QC fluid is not contaminated.</p> <p>Check the test setup to see if you need to change the parameters.</p> <p>Repeat the test in another way.</p> <p>Check the incubation temperature. Retest with a new reagent or QC.</p> |
| 13) The sample is contaminated. | <p>The sampler is not properly connected, resulting in leakage.</p> <p>Properly clean the sampling probe or replace the probe.</p> <p>Clean the wash tank.</p> <p>Ensure that there is no residual fluid on the probe.</p> <p>Check if the detergent is fresh and not dirty.</p> <p>When setting up the test, select the double rinse.</p> <p>Sample capacity is contaminated. Replace it with a new one.</p> |
| 14) The incubation cell is not hot. | <p>Check the temperature setting and set it to 37 °C.</p> <p>Check if the plug is plugged in.</p> |
| 15) Fluid leaks under the analyzer. | <p>The drain pipe was not inserted into the waste liquid tank.</p> |
| The analyzer does not turn on. | <p>The tank for waste liquid is full.</p> <p>Wash tank locked.</p> |

Warning: (1) In case of device failure, contact an Authorized Representative in the Russian Federation immediately to get technical support!

(2) Only specialists certified by the manufacturer can be allowed to repair the analyzer. To replace the spare parts, contact your dealer.

(3) When the analyzer is out of order due to repair or is intended for disposal, clean the fluid system, clean and sterilize the analyzer, remove the remaining fluid in different parts of the device with a dry cloth so that the instrument is dry. Ensure secure fixing of moving parts, such as a sampling probe, a cuvette grabber, etc.

(4) When repacking and transporting the analyzer, transportation should be carried out at a temperature of $-20^{\circ}\text{C} \sim 55^{\circ}\text{C}$, relative humidity $\leq 85\%$, barometric pressure $50\text{ kPa} \sim 106\text{ kPa}$. During transport, pay attention to the characters TOP, Fragile, afraid of dampness on the packaging and transport the device in accordance with the requirements.

(5) When repacking and transporting the device, take care to protect the front panel, the sampling probe, the cuvette grip unit, etc. as well as liquid connections on the rear panel from external force or contact with other objects or damage. During transport, the analyzer must remain upright and must not be tilted or placed on its side. Try to avoid vibration during operation. After service or repair, check it before use.

Annex I: Names and contents of toxic / hazardous substances or elements**1. Name and content of toxic / hazardous substances or elements**

Detail	Toxic / hazardous substances or elements					
	Lead (Pb)	Mercury (Hg)	Cadmium (CD)	Hexavalent chrome (Cr (VI))	Polybrominated biphenyls (PBB)	Polybrominated diphenyl ethers (PBDE)
Built-in circuit board	×	○	○	○	○	○
Housing	×	○	○	×	○	○
Display	×	○	○	○	○	○
Photovoltaic components	×	○	○	○	○	○
Internal electronic wire	○	○	○	○	○	○
Accessories	×	○	○	○	○	○
○: Indicates that the content of a toxic / hazardous substance in all homogeneous materials of a component is below the limit specified in SJ / T11363-2006. ×: Indicates the content of a toxic / hazardous substance in at least one homogeneous material of the component exceeds the limit specified in standard SJ / T11363-2006.						

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Sign value: An electronic information product contains certain toxic / hazardous substances that have a lifespan of 20 years. You can use it during this period. When this period expires, the product must enter the recycling system..

Chapter 16. Packing, storage, transportation

The device is packaged in a plywood box, ensuring the safety of the medical product from climatic and mechanical influences during transportation.

Inside the box, the device is surrounded on all sides by gaskets made of foamed polymer. Also, the foamed polymer gasket should be located inside the analyzer, protecting the compartments for reagents and samples, as well as the sampling and sampling cuvette units. The belts of the mechanisms for moving the sampler blocks and the cuvette gripping unit must be tightened with polymer clamps, which exclude accidental movement of the nodes.

Store the product in a dry place at a temperature from -20°C to $+55^{\circ}\text{C}$ and relative humidity $\leq 85\%$

Protect from moisture.

Avoid heat or storage under direct sunlight.

The research room should have an operating temperature of $+10^{\circ}\text{C}$ to $+30^{\circ}\text{C}$ and a relative humidity of $\leq 70\%$

If the appliance has cooled down during transport, please wait until it reaches room temperature before turning it on.

When putting into storage, the product must be packed in a supplier's (manufacturer's) packaging.

It is allowed to transport all types of transport in covered vehicles in accordance with the rules of transportation valid for this type of transport at a temperature range from -20°C to $+55^{\circ}\text{C}$ and relative air humidity up to 85%.

The device should be transported only in the packaging of the supplier (manufacturer).

The product may be transported in an additional shipping container of the transport company.

The method of packing shipping containers with products should exclude the possibility of their movement.

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In case of violation of the above conditions, the authorized representative in the Russian Federation of Laboratory Solutions LLC disclaims any responsibility for the quality of the results and for the operability of the device, and also has the right to refuse to replace the failed device or any part of it under warranty.

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