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Boschstraße 16  
D-21423 Winsen  
Germany

We, QualiMed Innovative Medizinprodukte GmbH, hereby confirm that the attached Instructions for Use for:

**Product Name**

- ATLAS Coronary Intravascular Stent Implantation System (Rapid Exchange)
- JUPITER Coronary Intravascular Stent Implantation System (Rapid Exchange)
- APOLLO Coronary Intravascular Stent Implantation System (Rapid Exchange)
- ABARIS Coronary Intravascular Stent Implantation System (Rapid Exchange)


**Accessories:**

- PEGASUS PTCA Balloon Catheter
- POLLUX Aspiration Catheter (Set)

are the original Instructions for Use.

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BIC-/SWIFT-Code: GENODE

## INSTRUCTIONS FOR USE

### ATLAS Coronary Intravascular Stent Implantation System (Rapid Exchange)

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout the instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

ATLAS Coronary Intravascular Stent Implantation System/ Rapid Exchange is a flexible balloon expandable stent mounted on a rapid exchange (RX) stent implantation system. The RX stent implantation system is designed for safe insertion of the stent in the coronary vessels and for stent deployment in the desired area by means of expansion of the balloon on which the stent is mounted.

The RX Stent Implantation System consists of a dilatation catheter with a balloon on the distal end. The stent is mounted on the balloon. The balloon of stent implantation system expands to the nominal diameter at the nominal pressure.

ATLAS Coronary Intravascular Stent Implantation System/ Rapid Exchange has the following attributes:

- Stent is made from cobalt-chromium alloy (L605).
- The balloon of the dilatation catheter of stent implantation system is made from nylon Grilamid L25.
- Balloon diameter – from 2 to 4 mm; balloon length from 12 to 40 mm.
- Tip of the dilatation catheter and flexible distal shaft are made from High-Density polyethylene (HDPE) Pebax 7233, and MEVOPUR Black colorant is added to the tip of the dilatation catheter.
- The dilatation catheter of stent implantation system is fitted with a three-layer inner tube made from High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant.
- In order to ensure best sliding and stent implantation to the target area, the distal shaft of the dilatation catheter of stent implantation system is fitted with two-layer hydrophilic coating. Base coat - DSM Primer 41002, Top coat - DSM Comfort Coat Topcoat 43003.
- Hypo tube of the dilatation catheter of stent implantation system consists of two parts: 1) (PTFE) Polytetrafluoroethylene (Xylan) coated stainless steel 304 L and 2) uncoated stainless steel 304 L.
- Kinking protection cover of the dilatation catheter of stent implantation system is made from Thermoplastic Urethane Texin RxT85A with Blue titanium TR-27 colorant.
- The cannula Luer Lock is made from Polycarbonate (PC) Lexan HPS 7. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint.
- The balloon of the dilatation catheter of stent implantation system has two radiopaque markers (Platinum – 90%, Iridium – 10%), used for precise positioning of the balloon.

- The balloon of stent implantation system can be used for post-dilatation purposes.
- The dilatation catheter of stent implantation system is fitted with two shaft markers made from stainless steel 304 L, designed to determine the position of the dilatation catheter implantation system in relation to the guiding catheter.
- Balloon nominal pressure – 9 bar.
- Burst pressure – 18 bar (for balloon diameter from 2 to 3.50 mm); 16 bar (for balloon diameter 4.00 mm).
- Dilatation catheter of stent implantation system is compatible with guide wire diameter 0.014”.
- The tip of dilatation catheter of stent implantation system is put on the protective guide wire made from stainless steel WKN 1.4301 (tip of protective guide wire from Polytetrafluoroethylene (PTFE) Teflon with Blue basic colorant), length 360, diameter 0.38 mm, passing on the inside of the distal part of the catheter and comes outside through the port.
- In order to protect the stent and the balloon from damage, on the side of the tip the dilatation catheter is fitted with a transparent cover (cap) made from Polytetrafluoroethylene (PTFE) Teflon, which should be removed before the stenting procedure along with retrievable protective guide wire.

**Table 1. Dimensions of ATLAS Coronary Intravascular Stent, dimensions of the balloon of dilatation catheter of stent implantation system**

| Stent diameter / Balloon diameter, mm | Stent length, mm | Balloon length, mm | Usable balloon length, $\pm 0.5$ mm | Inner distance of the markers, $\pm 0.5$ mm | Strut thickness, mm | Strut breadth, mm |
|---------------------------------------|------------------|--------------------|-------------------------------------|---------------------------------------------|---------------------|-------------------|
| 25                                    | 10               | 12                 | 12,50                               | 11,30                                       | 0,080 $\pm$ 0,010   | 0,104 $\pm$ 0,010 |
|                                       | 14               | 16                 | 16,50                               | 15,30                                       |                     |                   |
|                                       | 18               | 20                 | 20,50                               | 19,30                                       |                     |                   |
|                                       | 24               | 25                 | 25,50                               | 24,30                                       |                     |                   |
|                                       | 28               | 30                 | 30,50                               | 29,30                                       |                     |                   |
|                                       | 34               | 35                 | 35,50                               | 34,30                                       |                     |                   |
|                                       | 38               | 40                 | 41,50                               | 40,30                                       |                     |                   |
| 25                                    | 10               | 12                 | 12,50                               | 11,30                                       | 0,080 $\pm$ 0,010   | 0,104 $\pm$ 0,010 |
|                                       | 14               | 16                 | 16,50                               | 15,30                                       |                     |                   |
|                                       | 18               | 20                 | 20,50                               | 19,30                                       |                     |                   |
|                                       | 24               | 25                 | 25,50                               | 24,30                                       |                     |                   |
|                                       | 28               | 30                 | 30,50                               | 29,30                                       |                     |                   |
|                                       | 34               | 35                 | 35,50                               | 34,30                                       |                     |                   |
|                                       | 38               | 40                 | 41,50                               | 40,30                                       |                     |                   |
| 50                                    | 10               | 12                 | 12,50                               | 11,30                                       | 0,080 $\pm$ 0,010   | 0,104 $\pm$ 0,010 |
|                                       | 14               | 16                 | 16,50                               | 15,30                                       |                     |                   |
|                                       | 18               | 20                 | 20,50                               | 19,30                                       |                     |                   |
|                                       | 24               | 25                 | 25,50                               | 24,30                                       |                     |                   |

|   |    |    |       |       |               |               |
|---|----|----|-------|-------|---------------|---------------|
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
| 5 | 14 | 16 | 16,50 | 15,30 | 0,080 ± 0,010 | 0,104 ± 0,010 |
|   | 18 | 20 | 20,50 | 19,30 |               |               |
|   | 24 | 25 | 25,50 | 24,30 |               |               |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
|   | 14 | 16 | 16,50 | 15,30 |               |               |
| 0 | 18 | 20 | 20,50 | 19,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 24 | 25 | 25,50 | 24,30 |               |               |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
|   | 14 | 16 | 16,50 | 15,30 |               |               |
|   | 18 | 20 | 20,50 | 19,30 |               |               |
| 5 | 24 | 25 | 25,50 | 24,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 |               |               |
|   | 18 | 20 | 20,50 | 19,30 |               |               |
|   | 24 | 25 | 25,50 | 24,30 |               |               |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
| 0 | 14 | 16 | 16,50 | 15,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 18 | 20 | 20,50 | 19,30 |               |               |
|   | 24 | 25 | 25,50 | 24,30 |               |               |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
|   | 14 | 16 | 16,50 | 15,30 |               |               |
| 0 | 18 | 20 | 20,50 | 19,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 24 | 25 | 25,50 | 24,30 |               |               |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |
|   | 10 | 12 | 12,50 | 11,30 |               |               |
|   | 14 | 16 | 16,50 | 15,30 |               |               |
|   | 18 | 20 | 20,50 | 19,30 |               |               |
| 0 | 24 | 25 | 25,50 | 24,30 | 0,095 ± 0,010 | 0,099 ± 0,010 |
|   | 28 | 30 | 30,50 | 29,30 |               |               |
|   | 34 | 35 | 35,50 | 34,30 |               |               |
|   | 38 | 40 | 41,50 | 40,30 |               |               |

**Table 2. Technical characteristics of the dilatation catheter of stent implantation system**

| Parameter                                                                                  | Value                |
|--------------------------------------------------------------------------------------------|----------------------|
| Available length of dilatation catheter of stent implantation system, mm                   | 1415 ± 15            |
| Length of dilatation catheter of stent implantation system, mm                             | 3,5 ± 0,5            |
| Distal part of catheter, mm                                                                | 237±15               |
| Flexible distal shaft (Transition zone from hypo tube to outer tubing of the catheter), mm | 105±15               |
| Hypo tube, mm                                                                              | 1070±15              |
| Proximal shaft diameter of the dilatation catheter of stent                                | Max. 2,1 F (0,70 mm) |

|                                                               |                                                   |
|---------------------------------------------------------------|---------------------------------------------------|
| plantation system, F / mm                                     |                                                   |
| statal shaft diameter of the dilatation catheter of stent     | Max. 2,6 F (0,86 mm)                              |
| plantation system, F / mm                                     |                                                   |
| p - Proximal Shaft marker Distance of dilatation catheter, mm | 985 ± 15                                          |
| p - Distal Shaft marker Distance of dilatation catheter, mm   | 880 ± 15                                          |
| anking protection cover, mm                                   | Diameter – 1,10<br>Length - 35<br>Thickness – 0,2 |
| annula Luer Lock, mm                                          | Length - 41<br>Inner diameter from 0,7 to 1,25    |
| total balloon length of the dilatation catheter of stent      | Depends on stent size*                            |
| plantation system, mm                                         |                                                   |
| asable balloon length of the dilatation catheter of stent     | Depends on stent size*                            |
| plantation system, mm                                         |                                                   |
| ner distance of the markers                                   | Depends on stent size*                            |
| pered Tip Profile of the dilatation catheter of stent         | 0,45                                              |
| plantation system, mm                                         |                                                   |
| p entrance profile of the dilatation catheter of stent        | 0,54                                              |
| plantation system, mm                                         |                                                   |

#### ications for use

ATLAS Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for treating of coronary arteries when treating coronary vessel constriction due to advancing arterial sclerotic disease, not requiring open surgery.

ATLAS Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for improving the lumen of native coronary arteries (de novo) or vein grafts:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;
- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

The treated lesion length should be less than the nominal stent length (stent length - from 10 to 38 mm) with reference vessel diameters from 2.0 to 4.0 mm.

## **Contraindications**

Implantation of ATLAS intravascular coronary stent implantation system/rapid exchange is contraindicated for the patients having the following conditions:

- Allergic reaction to cobalt-chromium.
- Severe reaction to contrast agents.
- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.
- In-stent restenosis.
- Myocardial infarction < 72 hours.
- Stenting of saphenous vein grafts.
- Unprotected left main coronary artery.
- Total occlusion of target vessel.
- Heavily calcified lesions.
- Lesions involving arterial segments with highly tortuous anatomy.
- Lesions involving a bifurcation. Stent implantation in bifurcation area.
- Left ventricular ejection fraction < 30 %.
- Cardiogenic shock.
- Presence of definite or probable intraluminal thrombus.
- Any patients judged to have a lesion which may interfere with proper stent deployment.
- Pregnancy and breast-feeding.
- Patients with a history of gastrointestinal hemorrhage.
- Surgical treatment within two weeks before implantation, delivery, blood vessel biopsy and puncture of the vessel not susceptible to compression.
- Patients with recent acute cerebrovascular event (ACE).
- Diabetic hemorrhagic retinopathy and other conditions associated with coagulation defects.

## **Potential Adverse Events**

- Abrupt stent closure.
- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.
- Stenocardia.
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).
- Arteriovenous fistula.
- Cardiac tamponade.
- Cardiogenic shock.
- Death.

- Emboli, distal (air, tissue or thrombotic emboli).
- Emergent coronary artery bypass surgery.
- Heart failure.
- Hematoma.
- Haemorrhage.
- Hypotension/Hypertension.
- Infectious complications.
- Infectious inflammation.
- Myocardial ischemia.
- Coronary artery dissection, perforation or rupture.
- Pericardial effusion.
- Pseudo aneurysm, femoral.
- Renal failure.
- Restenosis of stented segment.
- Shock/Pulmonary edema.
- Coronary artery spasm.
- Stent embolization.
- Stent thrombosis/occlusion.
- Stroke/cerebrovascular accident/transient ischemic attack.
- Total occlusion of coronary artery or shunt.
- Vessel trauma requiring surgical repair or reintervention.
- Histologic changes in vessel wall, including inflammation, cellular damage or necrosis.
- Allergic reaction to cobalt-chromium.

## **arnings**

- Application of ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange carries an associated risk of subacute thrombosis, vascular complications such as bleeding, hematoma or pseudoaneurysm. During the entire vessel endothelialization period all patients should undergo anti-platelet therapy in the stenting area.
- ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange must not be implanted in patients with allergic reactions to cobalt-chromium.
- Implantation of ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange should be performed only by physicians who have received appropriate training.
- Implantation of ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.

- Direct Stenting without prior dilation of the lesion has not been studied using this product.
- Implantation of ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange may lead to development of subsequent restenosis, which may require reintervention (particularly, insertion and expansion of a dilatation catheter in arterial segment containing the stent). The long-term outcome following repeated dilation of coronary stents is unknown at present.
- When multiple stents are required, if placement results in stent to stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.

## Application method

### Precaution Prior to Use

Carefully inspect the sterile package before opening. Do not use ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange after the "Expire" date. If the integrity of the sterile package has been compromised prior to the "Expire" date (e.g., damage of the package), contact manufacturer QualiMed Innovative Medical Products GmbH, Germany or its authorized representative in the Russian Federation - CJSC "QualiMed". Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use. Do not use if any defects are noted in the device or package.

**NOTE:** If at any time during use the proximal shaft for catheter delivery has been bent or kinked, do not continue to use the device.

**Devices and materials required (not included in ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange package)** (All devices and materials are disposable: do not reuse or resterilize them):

- Guiding catheter(s) (inner lumen min. 1.42 mm) with the appropriate configuration for the coronary artery to be treated;
- 10 or 20 ml syringe;
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/ 2.39mm;
- Guide wire 0.014 inch/ 0.36 mm;
- 60% contrast medium diluted 1:1;
- Sterile physiological saline solution;
- Inflation device and a three-way stopcock for balloon inflation;
- Introducer.

## Preparation for Use

### Unboxing Removal (Step Action)

1. Carefully remove ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange from its protective tubing for preparation of the implantation system. Do not bend or kink hypo tube of the dilatation catheter of implantation system during removal.
2. Remove the transparent cover (cap) from the dilatation catheter from the tip and carefully retrieve the protective guide wire. Holding the dilatation catheter of implantation system with one hand on the proximal part from the stent, with the other hand carefully remove the transparent cover (cap) of the stent distally.

**NOTE:** Care should be taken not to kink or bend the shaft of the dilatation catheter upon application.

### **Stent System Guide wire Lumen Flush (Step Action)**

1. Flush stent system guide wire lumen with heparinized saline using a syringe.
2. Verify that the stent is positioned between the proximal and distal radiopaque balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted.

### **Preparation of the dilatation catheter of implantation system (Step Action)**

Prepare inflation device/syringe with diluted contrast medium.

Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypo tube when connecting to inflation device/syringe.

With tip down, orient stent implantation system vertically.

Open stopcock to stent implantation system; pull negative for 15 seconds; release to neutral for contrast fill.

Close stopcock to stent implantation system; purge inflation device/syringe of all air.

Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use device.

If a syringe was used, attach a prepared inflation device to stopcock.

Open stopcock to the dilatation catheter of stent implantation system.

Leave on neutral.

### **Stent Delivery Procedure (Step Action)**

1. Prepare the vascular access site according to standard percutaneous transluminal coronary angioplasty (PTCA) practice.
2. Predilate the lesion/vessel with appropriate diameter balloon having a ratio of 1:1 with the diameter of the vessel.
3. Maintain neutral pressure on inflation device attached to the dilatation catheter of stent implantation system.

4. Put the dilatation catheter of stent implantation system onto the proximal part of guide wire, while maintaining the distal part of guide wire position across target lesion.
5. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
6. Ensure that the guiding catheter is securely placed in the lesion area, then insert the dilatation catheter of implantation system in the coronary artery. Carefully advance the stent implantation system into the hub of the guiding catheter, keeping the hypo tube of the dilatation catheter of stent implantation system straight.

## TE

If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance indicate a problem and use of excessive force may result in stent damage or stent dislodgement from dilatation catheter of stent implantation system. Maintain guide wire placement across the lesion and move the stent implantation system and guiding catheter as a single unit.

7. Advance the stent implantation system over the guide wire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point for correct positioning. The inside edges of the marker bands indicate both the stent edges and balloon shoulders inflated. If the position of the stent is not optimal, do not deploy the stent; in such case the stent should be carefully repositioned or removed.
8. Sufficiently tighten the rotating hemostatic valve. Stent is now ready to be deployed.

## nt Deployment Procedure (Step Action)

1. Inflate the balloon of catheter of stent implantation system to a minimum pressure of 9 bar, (900 kPa) (deployment pressure). Higher pressures may be necessary to expand the stent to optimize stent apposition against the arterial wall. The pressure of the balloon of catheter of stent implantation system must not exceed rated burst pressure of 18 bar (1800 kPa) for stent diameters 2.00 to 3.50 mm and 16 bar (1600 kPa) for stent diameters at 4.00 mm.
  2. Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
  3. Deflate balloon by pulling negative on inflation device until balloon is fully deflated.
  4. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent.
- Optimal expansion requires that the stent be in full contact with the artery wall. All efforts should be taken to assure that the stent is not underdilated. Stent expansion should be performed with fluoroscopic visualization.

5. If stent sizing/apposition requires optimization, readvance the balloon of the catheter of stent implantation system, or another dilatation catheter of the appropriate size, to the stented area using standard angioplasty techniques.
6. Inflate the balloon to the desired pressure while observing by X-ray imaging. Deflate the balloon (see Balloon Compliance Chart supplied with device).
7. Reconfirm stent position and angiographic result. Ensure that the stent is positioned between the proximal and distal radiopaque markers of the balloon of stent implantation system. Repeat inflations until the desired result is achieved.

#### **Stent Implantation System Removal Procedure (Step Action)**

1. Ensure balloon of the catheter is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position and negative pressure on inflation device, withdraw implantation system.

#### **Precautions when handling ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Before using, check the product "Expire" date.
- Verify visually that the stent is located between the proximal and distal balloon markers of the implantation system.
- Coronary Intravascular Stent Implantation System/Rapid Exchange is designed for use as a unit. The stent is not to be removed from its delivery balloon, it may result in damaging and/or embolization of the stent. The stent is not designed to be crimped onto another balloon.
- Special care must be taken not to handle or in any way disrupt the stent position on the implantation system. This is most important during catheter removal from packaging, placement over guide wire, and advancement through hemostasis valve adapter and guiding catheter hub.
- Excessive manipulation, e.g. rolling the mounted stent, may result in dislodgement of the stent from the implantation system.
- Avoid handling the stent on the delivery hands without using protective gloves.
- In the event the stent is not deployed on the implantation system, please contact QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the RF CJSC "Lancet" to clarify the return procedure. Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use.
- Use only the appropriate solution of contrast medium. Do not use air or any gas medium to inflate the balloon.

Stent contact with any fluid prior to placement is not recommended. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

## **Precautions when implanting and removing ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange**

Do not pre-inflate balloon prior to stent deployment other than as directed.

Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).

The target lesion must be sufficiently predilated prior to stent implantation.

Do not expand the stent if it is not properly positioned in the vessel.

Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the Balloon Compliance Chart fixed to product label. Use of pressures higher than specified on product label may result in a ruptured balloon and potential intimal damage and dissection. The vessel should be pre-dilated with appropriate diameter balloon having a 1:1 ratio with the vessel diameter.

Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.

If unusual resistance is felt at any time during lesion access before stent implantation, the stent implantation system and the guiding catheter should be removed as a single unit.

Do not attempt to pull an unexpanded stent back into the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.

An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.

### **E: The following steps should be executed under direct visualization using X-ray imaging.**

Maintain guide wire placement across the lesion during the entire removal process. Carefully pull back the stent implantation system until the proximal radiopaque balloon marker of the stent implantation system is aligned with the distal tip of the guiding catheter.

The stent implantation system and guiding catheter should be pulled back until the tip of the guiding catheter is just distal to the arterial sheath, allowing the guiding catheter to straighten.

Carefully retract the stent implantation system into the guiding catheter and remove the stent implantation system and guiding catheter from the patient as a single unit while leaving the guide wire across the target lesion.

Failure to follow these steps and/or applying excessive force to the stent implantation system can potentially result in stent damage, stent dislodgement from the balloon and/or damage to the implantation system.

## **Precautions after implantation of ATLAS Coronary Intravascular Stent**

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, or a dilatation catheter to avoid disrupting the stent geometry.
- Do not perform Magnetic Resonance Imaging (MRI) scan on patient's post-stent implantation until the stent has been completely endothelialized (90 days) to minimize the potential for migration. The stent may cause artefacts in MRI scans due to distortion of the magnetic field.
- In order to reduce the risk of thrombosis inside the stent, during the entire vessel endothelialization period, anti-platelet therapy in the stenting area should be prescribed for all patients.

## **Single patient use precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

## **Single patient package content:**

1 ATLAS Coronary Intravascular Stent Implantation System/Rapid Exchange in dispenser with balloon compliance scale packed in sealed sterile plastic package. Package with device, instructions for use and desiccant placed in a cardboard box.

## **Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.
2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

## Disposal instructions

After use, dispose of product and packaging in accordance with common requirements to disposal of waste of medical devices.

## Shelf Life

Shelf life – 2 years from the date of sterilization. Do not use after expiry.


## Warranty/Liability

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the Instruction of Use. The warning notices being contained in this Instruction of Use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.


Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or the device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Postfachstrasse 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, Building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

## Graphical symbols for labeling primary and secondary packages

| Symbol                                                                              | Description                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Sterilized with ethylene oxide gas                                     |
|  | Lot number                                                             |
|  | Manufacturer                                                           |
|  | Each device is suitable for one (1) time use only, reuse is PROHIBITED |

|                                                                                                                                                                        |                                    |
|------------------------------------------------------------------------------------------------------------------------------------------------------------------------|------------------------------------|
|                                                                                        | Expires                            |
| <br> | Keep in a cool, dry and dark place |
|                                                                                       | Read instructions prior to use     |
|                                                                                       | Hydrophilic coating                |
|                                                                                     | CE Marking                         |
|                                                                                     | Logo:<br>IK /INTELLIGENT KINETICS  |

**Manufacturer**  
**aliMed Innovative Medizinprodukte GmbH**  
 Schwanenstraße 16  
 423 Witten  
 Germany  
 Phone No.: +49 4171 6578 0  
 Fax: +49 4171 6578 11

## INSTRUCTIONS FOR USE

### JUPITER Coronary Intravascular Stent Implantation System (Rapid Exchange)

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout the instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

JUPITER Coronary Intravascular Stent Implantation System/ Rapid Exchange is a flexible balloon expandable stent mounted on a rapid exchange (RX) stent implantation system. The RX stent implantation system is designed for safe insertion of the stent in the coronary vessels and for stent deployment in the desired area by means of expansion of the balloon on which the stent is mounted.

The RX Stent Implantation System consists of a dilatation catheter with a balloon on the distal end. The stent is mounted on the balloon. The balloon of stent implantation system expands to the nominal diameter at the nominal pressure.

JUPITER Coronary Intravascular Stent Implantation System/ Rapid Exchange has the following features:

- Stent is made from stainless steel 316 LVM.
- The balloon of the dilatation catheter of stent implantation system is made from nylon Grilamid L25.
- Balloon diameter – from 2 to 4 mm; balloon length from 11 to 40 mm.
- Tip of the dilatation catheter and flexible distal shaft are made from High-Density polyethylene (HDPE) Pebax 7233, and MEVOPUR Black colorant is added to the tip of the dilatation catheter.
- The dilatation catheter of stent implantation system is fitted with a three-layer inner tube made from High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant.
- In order to ensure best sliding and stent implantation to the target area, the distal shaft of the dilatation catheter of stent implantation system is fitted with two-layer hydrophilic coating. Base coat - DSM Primer 41002, Top coat - DSM Comfort Coat Topcoat 43003.
- Hypo tube of the dilatation catheter of stent implantation system consists of two parts: 1) (PTFE) Polytetrafluoroethylene (Xylan) coated stainless steel 304 L and 2) uncoated stainless steel 304 L.
- Kinking protection cover of the dilatation catheter of stent implantation system is made from Thermoplastic Urethane Texin RxT85A with Blue titanium TR-27 colorant.
- The cannula Luer Lock is made from Polycarbonate (PC) Lexan HPS 7. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint.

- The balloon of the dilatation catheter of stent implantation system has two radiopaque markers (Platinum – 90%, Iridium– 10%), used for precise positioning of the balloon.
- The balloon of stent implantation system can be used for post-dilatation purposes.
- The dilatation catheter of stent implantation system is fitted with two shaft markers made from stainless steel 304 L, designed to determine the position of the dilatation catheter implantation system in relation to the guiding catheter.
- Balloon nominal pressure – 9 bar.
- Burst pressure – 18 bar (for balloon diameter from 2 to 3.50 mm); 16 bar (for balloon diameter 4.00 mm).
- Dilatation catheter of stent implantation system is compatible with guide wire diameter 0.014”.
- The tip of dilatation catheter of stent implantation system is put on the protective guide wire made from stainless steel WKN 1.4301 (tip of protective guide wire from Polytetrafluoroethylene (PTFE) Teflon with Blue basic colorant), length 360, diameter 0.38 mm, passing on the inside of the distal part of the catheter and comes outside through the port.
- In order to protect the stent and the balloon from damage, on the side of the tip the dilatation catheter is fitted with a transparent cover (cap) made from Polytetrafluoroethylene (PTFE) Teflon, which should be removed before the stenting procedure along with retrievable protective guide wire.

**Table 1. Dimensions of JUPITER Coronary Intravascular Stent, dimensions of the balloon of dilatation catheter of stent implantation system**

| Stent diameter / Balloon diameter, mm | Stent length, mm | Balloon length, mm | Usable balloon length, ± 0.5 mm | Inner distance of the markers, ± 0.5 mm | Strut thickness, mm | Strut breadth, mm |
|---------------------------------------|------------------|--------------------|---------------------------------|-----------------------------------------|---------------------|-------------------|
| 25                                    | 10               | 11                 | 11,50                           | 10,3                                    | 0,105 ± 0,010       | 0,118 ± 0,010     |
|                                       | 14               | 15                 | 15,50                           | 14,3                                    |                     |                   |
|                                       | 18               | 20                 | 20,50                           | 19,3                                    |                     |                   |
|                                       | 24               | 25                 | 25,50                           | 24,3                                    |                     |                   |
|                                       | 28               | 30                 | 30,50                           | 29,3                                    |                     |                   |
|                                       | 34               | 35                 | 35,50                           | 34,3                                    |                     |                   |
|                                       | 38               | 40                 | 41,50                           | 40,3                                    |                     |                   |
| 25                                    | 10               | 11                 | 11,50                           | 10,3                                    | 0,105 ± 0,010       | 0,118 ± 0,010     |
|                                       | 14               | 15                 | 15,50                           | 14,3                                    |                     |                   |
|                                       | 18               | 20                 | 20,50                           | 19,3                                    |                     |                   |
|                                       | 24               | 25                 | 25,50                           | 24,3                                    |                     |                   |
|                                       | 28               | 30                 | 30,50                           | 29,3                                    |                     |                   |
|                                       | 34               | 35                 | 35,50                           | 34,3                                    |                     |                   |
|                                       | 38               | 40                 | 41,50                           | 40,3                                    |                     |                   |
| 50                                    | 10               | 11                 | 11,50                           | 10,3                                    | 0,105 ± 0,010       | 0,118 ± 0,010     |
|                                       | 14               | 15                 | 15,50                           | 14,3                                    |                     |                   |

|     |    |    |       |      |               |               |
|-----|----|----|-------|------|---------------|---------------|
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |
| 5   | 10 | 11 | 11,50 | 10,3 | 0,105 ± 0,010 | 0,118 ± 0,010 |
|     | 14 | 15 | 15,50 | 14,3 |               |               |
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |
| 10  | 10 | 11 | 11,50 | 10,3 | 0,105 ± 0,010 | 0,118 ± 0,010 |
|     | 14 | 15 | 15,50 | 14,3 |               |               |
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |
| 25  | 10 | 12 | 12,50 | 11,3 | 0,115 ± 0,010 | 0,125 ± 0,010 |
|     | 14 | 16 | 16,50 | 15,3 |               |               |
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |
| 50  | 10 | 12 | 12,50 | 11,3 | 0,115 ± 0,010 | 0,125 ± 0,010 |
|     | 14 | 16 | 16,50 | 15,3 |               |               |
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |
| 100 | 10 | 12 | 12,50 | 11,3 | 0,115 ± 0,010 | 0,125 ± 0,010 |
|     | 14 | 16 | 16,50 | 15,3 |               |               |
|     | 18 | 20 | 20,50 | 19,3 |               |               |
|     | 24 | 25 | 25,50 | 24,3 |               |               |
|     | 28 | 30 | 30,50 | 29,3 |               |               |
|     | 34 | 35 | 35,50 | 34,3 |               |               |
|     | 38 | 40 | 41,50 | 40,3 |               |               |

**Table 2. Technical characteristics of the dilatation catheter of stent implantation system**

| Parameter                                                                                  | Value     |
|--------------------------------------------------------------------------------------------|-----------|
| Available length of dilatation catheter of stent implantation system, mm                   | 1415 ± 15 |
| Length of dilatation catheter of stent implantation system, mm                             | 3,5 ± 0,5 |
| Distal part of catheter, mm                                                                | 237±15    |
| Flexible distal shaft (Transition zone from hypo tube to inner tubing of the catheter), mm | 105±15    |
| Hypo tube, mm                                                                              | 1070±15   |

|                                                                                         |                                                   |
|-----------------------------------------------------------------------------------------|---------------------------------------------------|
| Proximal shaft diameter of the dilatation catheter of stent implantation system, F / mm | Max. 2,1 F (0,70 mm)                              |
| Distal shaft diameter of the dilatation catheter of stent implantation system, F / mm   | Max. 2,6 F (0,86 mm)                              |
| Proximal Shaft marker Distance of dilatation catheter, mm                               | 985 ± 15                                          |
| Distal Shaft marker Distance of dilatation catheter, mm                                 | 880 ± 15                                          |
| Stenting protection cover, mm                                                           | Diameter – 1,10<br>Length - 35<br>Thickness – 0,2 |
| Annula Luer Lock, mm                                                                    | Length - 41<br>Inner diameter from 0,7 to 1,25    |
| Total balloon length of the dilatation catheter of stent implantation system, mm        | Depends on stent size*                            |
| Deployable balloon length of the dilatation catheter of stent implantation system, mm   | Depends on stent size*                            |
| Marker distance of the markers                                                          | Depends on stent size*                            |
| Preperforated Tip Profile of the dilatation catheter of stent implantation system, mm   | 0,45                                              |
| Entrance profile of the dilatation catheter of stent implantation system, mm            | 0,54                                              |

#### Indications for use

JUPITER Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for treatment of coronary arteries when treating coronary vessel constriction due to advancing arterial sclerotic disease, not requiring open surgery.

JUPITER Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for improving the lumen of native coronary arteries (de novo) or vein grafts:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;
- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

The treated lesion length should be less than the nominal stent length (stent length - from 10 to 38 mm) with reference vessel diameters from 2.0 to 4.0 mm.

## **Contraindications**

Implantation of JUPITER intravascular coronary stent implantation system/rapid exchange is contraindicated for the patients having the following conditions:

- Allergic reaction to stainless steel.
- Allergic reaction to contrast agents.
- Patients with contraindications to antiplatelet and/or anticoagulant therapy.
- In-stent restenosis.
- Myocardial infarction < 72 hours.
- Stenting of saphenous vein grafts.
- Unprotected left main coronary artery.
- Total occlusion of target vessel.
- Heavily calcified lesions.
- Lesions involving arterial segments with highly tortuous anatomy.
- Lesions involving a bifurcation.
- Stent implantation in bifurcation area.
- Left ventricular ejection fraction < 30 %.
- Cardiogenic shock.
- Presence of definite or probable intraluminal thrombus.
- Any patients judged to have a lesion which may interfere with proper stent deployment.
- Pregnancy and breast-feeding.
- Patients with a history of gastrointestinal hemorrhage.
- Surgical treatment within two weeks before implantation, delivery, blood vessel biopsy and puncture of the vessel not susceptible to compression.
- Patients with recent acute cerebrovascular event (ACE).
- Diabetic hemorrhagic retinopathy and other conditions associated with coagulation defects.

## **Potential Adverse Events**

- Sudden stent closure.
- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.
- Angina pectoris.
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).
- Arteriovenous fistula.
- Cardiac tamponade.

- Cardiogenic shock.
- Death.
- Emboli, distal (air, tissue or thrombotic emboli).
- Emergent coronary artery bypass surgery.
- Heart failure.
- Hematoma.
- Haemorrhage.
- Hypotension/Hypertension.
- Infectious complications.
- Infectious inflammation.
- Myocardial ischemia.
- Coronary artery dissection, perforation or rupture.
- Pericardial effusion.
- Pseudo aneurysm, femoral artery.
- Renal failure.
- Restenosis of stented segment.
- Shock/Pulmonary edema.
- Coronary artery spasm.
- Stent embolization.
- Stent thrombosis/occlusion.
- Stroke/cerebrovascular accident/transient ischemic attack.
- Total occlusion of coronary artery or shunt.
- Vessel trauma requiring surgical repair or reintervention.
- Histological changes in vessel wall, including inflammation, cellular damage or necrosis.
- Allergic reaction to stainless steel.

#### **Warnings**

- Application of JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange carries an associated risk of subacute thrombosis, vascular complications such as bleeding, hematoma or pseudoaneurysm. During the entire vessel endothelialization period all patients should undergo anti-platelet therapy in the stenting area.
- JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange must not be implanted in patients with allergic reactions to stainless steel.
- Implantation of JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange should be performed only by physicians who have received appropriate training.

- Implantation of JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- Direct Stenting without prior dilation of the lesion has not been studied using this product.
- Implantation of JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange may lead to development of subsequent restenosis, which may require reintervention (particularly, insertion and expansion of a dilatation catheter in arterial segment containing the stent). The long-term outcome following repeated dilation of coronary stents is unknown at present.
- When multiple stents are required, if placement results in stent to stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.

## Application method

### Inspection Prior to Use

Carefully inspect the sterile package before opening. Do not use JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange after the "Expire" date. If the integrity of the sterile package has been compromised prior to the product "Expire" date (e.g., damage of the package), contact manufacturer QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the Russian Federation - CJSC "Medica". Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use. Do not use if any defects are noted in the device or package.

**NOTE:** If at any time during use the proximal shaft for catheter delivery has been bent or kinked, do not continue to use the device.

**Devices and materials required (not included in JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange package)** (All devices and materials are disposable: do not reuse or resterilize them):

- Guiding catheter(s) (inner lumen min. 1.42 mm) with the appropriate configuration for the coronary artery to be treated;
- 10 or 20 ml syringe;
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/ 2.39mm;
- Guide wire 0.014 inch/ 0.36 mm;
- 60% contrast medium diluted 1:1;
- Sterile physiological saline solution;
- Inflation device and a three-way stopcock for balloon inflation;
- Introducer.

## Preparation for Use

### **Stent Packaging Removal (Step Action)**

1. Carefully remove JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange from its protective tubing for preparation of the implantation system. Do not bend or kink hypo tube of the dilatation catheter of implantation system during removal.
2. Remove the transparent cover (cap) from the dilatation catheter from the tip and carefully retrieve the protective guide wire. Holding the dilatation catheter of implantation system with one hand on the proximal part from the stent, with the other hand carefully remove the transparent cover (cap) of the stent distally.

**NOTE:** Care should be taken not to kink or bend the shaft of the dilatation catheter upon application.

### **Stent System Guide wire Lumen Flush (Step Action)**

1. Flush stent system guide wire lumen with heparinized saline using a syringe.
2. Verify that the stent is positioned between the proximal and distal radiopaque balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted.

### **Preparation of the dilatation catheter of implantation system (Step Action)**

Prepare inflation device/syringe with diluted contrast medium.

Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypo tube when connecting to inflation device/syringe.

With tip down, orient stent implantation system vertically.

Open stopcock to stent implantation system; pull negative for 15 seconds; release to neutral for contrast fill.

Close stopcock to stent implantation system; purge inflation device/syringe of all air.

Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use device.

If a syringe was used, attach a prepared inflation device to stopcock.

Open stopcock to the dilatation catheter of stent implantation system.

Leave on neutral.

### **Stent Delivery Procedure (Step Action)**

1. Prepare the vascular access site according to standard percutaneous transluminal coronary angioplasty (PTCA) practice.
2. Predilate the lesion/vessel with appropriate diameter balloon having a ratio of 1:1 with the diameter of the vessel.

5. Maintain neutral pressure on inflation device attached to the dilatation catheter of stent implantation system.
6. Put the dilatation catheter of stent implantation system onto the proximal part of guide wire, while maintaining the distal part of guide wire position across target lesion.
7. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
8. Ensure that the guiding catheter is securely placed in the lesion area, then insert the dilatation catheter of implantation system in the coronary artery. Carefully advance the stent implantation system into the hub of the guiding catheter, keeping the hypo tube of the dilatation catheter of stent implantation system straight.

TE

If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance indicate a problem and use of excessive force may result in stent damage or stent dislodgement from dilatation catheter of stent implantation system. Maintain guide wire placement across the lesion and move the stent implantation system and guiding catheter as a single unit.

7. Advance the stent implantation system over the guide wire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point for correct positioning. The inside edges of the marker bands indicate both the stent edges and balloon shoulders inflated. If the position of the stent is not optimal, do not deploy the stent; in such case the stent should be carefully repositioned or removed.
8. Sufficiently tighten the rotating hemostatic valve. Stent is now ready to be deployed.

### Stent Deployment Procedure (Step Action)

1. Inflate the balloon of catheter of stent implantation system to a minimum pressure of 9 bar, (900 kPa) (deployment pressure). Higher pressures may be necessary to expand the stent to optimize stent apposition against the arterial wall. The pressure of the balloon of catheter of stent implantation system must not exceed rated burst pressure of 18 bar (1800 kPa) for stent diameters 2.00 to 3.50 mm and 16 bar (1600 kPa) for stent diameters at 4.00 mm.
2. Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
3. Deflate balloon by pulling negative on inflation device until balloon is fully deflated.
4. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent.

Optimal expansion requires that the stent be in full contact with the artery wall. All efforts should be taken to assure that the stent is not underdilated. Stent expansion should be performed with fluoroscopic visualization.

5. If stent sizing/apposition requires optimization, readvance the balloon of the catheter of stent implantation system, or another dilatation catheter of the appropriate size, to the stented area using standard angioplasty techniques.
6. Inflate the balloon to the desired pressure while observing by X-ray imaging. Deflate the balloon (see Balloon Compliance Chart supplied with device).
7. Reconfirm stent position and angiographic result. Ensure that the stent is positioned between the proximal and distal radiopaque markers of the balloon of stent implantation system. Repeat inflations until the desired result is achieved.

#### **Stent Implantation System Removal Procedure (Step Action)**

1. Ensure balloon of the catheter is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position and negative pressure on inflation device, withdraw implantation system.

#### **Precautions when handling JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Before using, check the product “Expire” date.
- Verify visually that the stent is located between the proximal and distal balloon markers of the implantation system.
- Coronary Intravascular Stent Implantation System/Rapid Exchange is designed for use as a unit. The stent is not to be removed from its delivery balloon, it may result in damaging and/or embolization of the stent. The stent is not designed to be crimped onto another balloon.
- Special care must be taken not to handle or in any way disrupt the stent position on the implantation system. This is most important during catheter removal from packaging, placement over guide wire, and advancement through hemostasis valve adapter and guiding catheter hub.
- Excessive manipulation, e.g. rolling the mounted stent, may result in dislodgement of the stent from the implantation system.
- Avoid handling the stent on the delivery hands without using protective gloves.
- In the event the stent is not deployed on the implantation system, please contact QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the RF CJSC

"Lancet" to clarify the return procedure. Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use.

- Use only the appropriate solution of contrast medium. Do not use air or any gas medium to inflate the balloon.
- Stent contact with any fluid prior to placement is not recommended. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

## **Precautions when implanting and removing JUPITER Coronary Intravascular Stent plantation System/Rapid Exchange**

- Do not pre-inflate balloon prior to stent deployment other than as directed.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- The target lesion must be sufficiently predilated prior to stent implantation.
- Do not expand the stent if it is not properly positioned in the vessel.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the Balloon Compliance Chart fixed to product label. Use of pressures higher than specified on product label may result in a ruptured balloon and potential intimal damage and dissection. The vessel should be pre-dilated with appropriate diameter balloon having a 1:1 ratio with the vessel diameter.
- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.
- If unusual resistance is felt at any time during lesion access before stent implantation, the stent implantation system and the guiding catheter should be removed as a single unit.
- Do not attempt to pull an unexpanded stent back into the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.

**NOTE:** The following steps should be executed under direct visualization using X-ray imaging.

- Maintain guide wire placement across the lesion during the entire removal process. Carefully pull back the stent implantation system until the proximal radiopaque balloon marker of the stent implantation system is aligned with the distal tip of the guiding catheter.
- The stent implantation system and guiding catheter should be pulled back until the tip of the guiding catheter is just distal to the arterial sheath, allowing the guiding catheter to straighten. Carefully retract the stent implantation system into the guiding catheter and remove the stent implantation system and guiding catheter from the patient as a single unit while leaving the guide wire across the target lesion.

Failure to follow these steps and/or applying excessive force to the stent implantation system can potentially result in stent damage, stent dislodgement from the balloon and/or damage to the implantation system.

### **Precautions after implantation of JUPITER Coronary Intravascular Stent**

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, or a dilatation catheter to avoid disrupting the stent geometry.
- Do not perform Magnetic Resonance Imaging (MRI) scan on patient's post-stent implantation until the stent has been completely endothelialized (90 days) to minimize the potential for migration. The stent may cause artefacts in MRI scans due to distortion of the magnetic field.
- In order to reduce the risk of thrombosis inside the stent, during the entire vessel endothelialization period, anti-platelet therapy in the stenting area should be prescribed for all patients.

### **Use precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Package content:** 1 JUPITER Coronary Intravascular Stent Implantation System/Rapid Exchange in dispenser with balloon compliance scale packed in sealed sterile plastic package. Package with device, instructions for use and desiccant placed in a cardboard box.

### **Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.

2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

### Disposal instructions

After use, dispose of product and packaging in accordance with common requirements to disposal waste of medical devices.

### Shelf Life

Shelf life – 2 years from the date of sterilization. Do not use after expiry.


### Warranty/Liability

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the Instruction of Use. The warning notices being contained in this Instruction of Use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.



Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Schwanenstrasse 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, Building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

### Graphical symbols for labeling primary and secondary packages

| Symbol                                                                              | Description                        |
|-------------------------------------------------------------------------------------|------------------------------------|
|  | Sterilized with ethylene oxide gas |
|  | Lot number                         |
|  | Manufacturer                       |

|                                                                                     |                                                                        |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|    | Each device is suitable for one (1) time use only, reuse is PROHIBITED |
|    | Expires                                                                |
|    | Keep in a cool, dry and dark place                                     |
|    | Read instructions prior to use                                         |
|  | Hydrophilic coating                                                    |
|  | CE Marking                                                             |
|  | Logo:<br>IK /INTELLIGENT KINETICS                                      |

**Manufacturer**  
**AliMed Innovative Medizinprodukte GmbH**  
 Schwanenbachstrasse 16  
 30523 Winsen  
 Lower Saxony, Germany  
 Phone No.: +49 4171 6578 0  
 Fax No.: +49 4171 6578 11

## INSTRUCTIONS FOR USE

### APOLLO Coronary Intravascular Stent Implantation System (Rapid Exchange)

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

APOLLO Coronary Intravascular Stent Implantation System/ Rapid Exchange is a flexible balloon expandable stent mounted on a rapid exchange (RX) stent implantation system. The RX stent implantation system is designed for safe insertion of the stent in the coronary vessels and for stent deployment in the desired area by means of expansion of the balloon on which the stent is mounted.

The RX Stent Implantation System consists of a dilatation catheter with a balloon on the distal end. The stent is mounted on the balloon. The balloon of stent implantation system expands to the nominal diameter at the nominal pressure.

APOLLO Coronary Intravascular Stent Implantation System/ Rapid Exchange has the following attributes:

- Stent is made from stainless steel (316 LVM) with biostable polymeric Paclitaxel-eluting coating. Depending on the stent size Paclitaxel content varies from 261 to 1289 µg. Paclitaxel concentration – 1.0 µg/mm<sup>2</sup>. Polymeric coating - ProTector (INN - Polysulfone).
- The balloon of the dilatation catheter of stent implantation system is made from nylon Grilamid L25.
- Balloon diameter – from 2 to 4 mm; balloon length from 11 to 40 mm.
- Tip of the dilatation catheter and flexible distal shaft are made from High-Density polyethylene (HDPE) Pebax 7233, and MEVOPUR Black colorant is added to the tip of the dilatation catheter.
- The dilatation catheter of stent implantation system is fitted with a three-layer inner tube made from High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant.
- In order to ensure best sliding and stent implantation to the target area, the distal shaft of the dilatation catheter of stent implantation system is fitted with two-layer hydrophilic coating. Base coat - DSM Primer 41002, Top coat - DSM Comfort Coat Topcoat 43003.
- Hypo tube of the dilatation catheter of stent implantation system consists of two parts: 1) (PTFE) Polytetrafluoroethylene (Xylan) coated stainless steel 304 L and 2) uncoated stainless steel 304 L.
- Kinking protection cover of the dilatation catheter of stent implantation system is made from Thermoplastic Urethane Texin RxT85A with Blue titanium dioxide TR-27 colorant.

- The cannula Luer Lock is made from Polycarbonate (PC) Lexan HPS 7. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint.
- The balloon of the dilatation catheter of stent implantation system has two radiopaque markers (Platinum – 90%, Iridium– 10%), used for precise positioning of the balloon.
- The balloon of stent implantation system can be used for post-dilatation purposes.
- The dilatation catheter of stent implantation system is fitted with two shaft markers made from stainless steel 304 L, designed to determine the position of the dilatation catheter implantation system in relation to the guiding catheter.
- Balloon nominal pressure – 9 bar.
- Burst pressure – 18 bar (for balloon diameter from 2 to 3.50 mm); 16 bar (for balloon diameter 4.00 mm).
- Dilatation catheter of stent implantation system is compatible with guide wire diameter 0.014”.
- The tip of dilatation catheter of stent implantation system is put on the protective guide wire made from stainless steel WKN 1.4301 (tip of protective guide wire from Polytetrafluoroethylene (PTFE) Teflon with Blue basic colorant), length 360, diameter 0.38 mm, passing on the inside of the distal part of the catheter and comes outside through the port.
- In order to protect the stent and the balloon from damage, on the side of the tip the dilatation catheter is fitted with a transparent cover (cap) made from Polytetrafluoroethylene (PTFE) Teflon, which should be removed before the stenting procedure along with retrievable protective guide wire.

**Table 1. Dimensions of APOLLO Coronary Intravascular Stent, dimensions of the balloon of the dilatation catheter of stent implantation system**

| Stent length, mm | Balloon length, mm | Usable balloon length, $\pm 0.5$ mm | Inner distance of the markers, $\pm 0.5$ mm | Strut thickness, mm | Strut breadth, mm |
|------------------|--------------------|-------------------------------------|---------------------------------------------|---------------------|-------------------|
| 10               | 11                 | 11,50                               | 10,30                                       | 0,105 $\pm$ 0,010   | 0,118 $\pm$ 0,010 |
| 14               | 15                 | 15,50                               | 14,30                                       |                     |                   |
| 18               | 20                 | 20,50                               | 19,30                                       |                     |                   |
| 24               | 25                 | 25,50                               | 24,30                                       |                     |                   |
| 28               | 30                 | 30,50                               | 29,30                                       |                     |                   |
| 34               | 35                 | 35,50                               | 34,30                                       |                     |                   |
| 38               | 40                 | 41,50                               | 40,30                                       |                     |                   |
| 10               | 11                 | 11,50                               | 10,30                                       | 0,105 $\pm$ 0,010   | 0,118 $\pm$ 0,010 |
| 14               | 15                 | 15,50                               | 14,30                                       |                     |                   |
| 18               | 20                 | 20,50                               | 19,30                                       |                     |                   |
| 24               | 25                 | 25,50                               | 24,30                                       |                     |                   |
| 28               | 30                 | 30,50                               | 29,30                                       |                     |                   |
| 34               | 35                 | 35,50                               | 34,30                                       |                     |                   |
| 38               | 40                 | 41,50                               | 40,30                                       |                     |                   |

|   |    |    |       |       |               |                |
|---|----|----|-------|-------|---------------|----------------|
| 0 | 10 | 11 | 11,50 | 10,30 | 0,105 ± 0,010 | 0,118 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 5 | 10 | 11 | 11,50 | 10,30 | 0,105 ± 0,010 | 0,118 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 11 | 11,50 | 10,30 | 0,105 ± 0,010 | 0,118 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 5 | 10 | 12 | 12,50 | 11,30 | 0,115 ± 0,010 | 0, 125 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,115 ± 0,010 | 0, 125 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,115 ± 0,010 | 0, 125 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |

**Table 2. Technical characteristics of the dilatation catheter of stent implantation system**

| Parameter                                                                        | Value     |
|----------------------------------------------------------------------------------|-----------|
| Available length of dilatation catheter of stent implantation system, mm         | 1415 ± 15 |
| Length of dilatation catheter of stent implantation system, mm                   | 3,5 ± 0,5 |
| Length of distal part of catheter, mm                                            | 237±15    |
| Length of flexible distal shaft (Transition zone from hypo tube to catheter), mm | 105±15    |

|                                                                                         |                                                   |
|-----------------------------------------------------------------------------------------|---------------------------------------------------|
| Proximal shaft diameter of the dilatation catheter of stent implantation system, F / mm | 1070±15                                           |
| Distal shaft diameter of the dilatation catheter of stent implantation system, F / mm   | Max. 2,1 F (0,70 mm)                              |
| Proximal Shaft marker Distance of dilatation catheter, mm                               | Max. 2,6 F (0,86 mm)                              |
| Distal Shaft marker Distance of dilatation catheter, mm                                 | 985 ± 15                                          |
| Working protection cover, mm                                                            | 880 ± 15                                          |
| Annula Luer Lock, mm                                                                    | Diameter – 1,10<br>Length - 35<br>Thickness – 0,2 |
| Total balloon length of the dilatation catheter of stent implantation system, mm        | Length - 41<br>Inner diameter from 0,7 to 1,25    |
| Expandable balloon length of the dilatation catheter of stent implantation system, mm   | Depends on stent size*                            |
| Marker distance of the markers                                                          | Depends on stent size*                            |
| Proximal Tip Profile of the dilatation catheter of stent implantation system, mm        | Depends on stent size*                            |
| Distal entrance profile of the dilatation catheter of stent implantation system, mm     | 0,45                                              |
|                                                                                         | 0,54                                              |

## Indications for use

APOLLO Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for treatment of coronary arteries when treating coronary vessel constriction due to advancing arterial sclerotic disease, not requiring open surgery.

APOLLO Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for improving the lumen of native coronary arteries (de novo) or vein grafts:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;
- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

Drug coating of APOLLO coronary intravascular stent significantly helps reduce:

- restenosis frequency after percutaneous coronary intervention,
- revascularization frequency of target lesion,
- frequency of late loss of vessel lumen after stenting procedure.

Implantation of APOLLO Drug-Eluting Coronary Intravascular Stent is also indicated for treatment of abrupt or threatening closure in patients with failed interventional therapy.

The treated lesion length should be less than the nominal stent length (stent length - from 10 to 38 mm) with reference vessel diameters from 2.0 to 4.0 mm.

### **Contraindications**

Implantation of APOLLO Coronary Intravascular Stent Implantation System/ Rapid Exchange is contraindicated for patients with the following conditions:

- Hypersensitivity to Paclitaxel.
- Allergic reaction to stainless steel.
- Allergic reaction to Polysulfone.
- Severe reaction to contrast agents.
- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.
- In-stent restenosis.
- Myocardial infarction < 72 hours.
- Stenting of saphenous vein grafts.
- Unprotected left main coronary artery.
- Total occlusion of target vessel.
- Heavily calcified lesions.
- Lesions involving arterial segments with highly tortuous anatomy.
- Lesions involving a bifurcation. Stent implantation in bifurcation area.
- Left ventricular ejection fraction < 30 %.
- Cardiogenic shock.
- Presence of definite or probable intraluminal thrombus.
- Any patients judged to have a lesion which may interfere with proper stent deployment.
- Pregnancy and breast-feeding.

### **Potential Adverse Events**

- Abrupt stent closure.
- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.

Stenocardia.

Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).

Arteriovenous Fistula.

Cardiac Tamponade.

Cardiogenic Shock.

Death.

Emboli, distal (air, tissue or thrombotic emboli).

Emergent Coronary Artery Bypass Surgery.

Heart Failure.

Hematoma.

Haemorrhage.

Hypotension/Hypertension.

Infectious complications.

Infectious inflammation.

Myocardial ischemia.

Coronary artery dissection, perforation or rupture.

Pericardial effusion.

Pseudo aneurysm, femoral.

Renal failure.

Restenosis of stented segment.

Shock/Pulmonary edema.

Coronary artery spasm.

Stent embolization.

Stent thrombosis/occlusion.

Stroke/cerebrovascular accident/transient ischemic attack.

Total occlusion of coronary artery or shunt.

Vessel trauma requiring surgical repair or reintervention.

Histologic changes in vessel wall, including inflammation, cellular damage or necrosis.

Allergic reaction to stainless steel.

Allergic reaction to Paclitaxel, Polysulfone.

Alopecia

Myalgia/Arthralgia

Gastrointestinal symptoms

Peripheral neuropathy

- Hematologic dyscrasia (anaemia, leukopenia, neutropenia, thrombocytopenia)
- Hepatic enzyme changes

## Warnings

- Application of APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange carries an associated risk of subacute thrombosis, vascular complications such as bleeding, hematoma or pseudoaneurysm. Antiplatelet therapy must be prescribed for all patients for a period of 6 months post procedure.
- APOLLO Stent Implantation System/Rapid Exchange must not be implanted in patients with allergic reactions to stainless steel, Polysulfone and/or hypersensitivity to Paclitaxel.
- Implantation of APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange should be performed only by physicians who have received appropriate training.
- Implantation of APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- Direct Stenting without prior dilation of the lesion has not been studied using this product.
- Implantation of APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange may lead to development of subsequent restenosis, which may require reintervention (particularly, insertion and expansion of a dilatation catheter in arterial segment containing the stent). The long-term outcome following repeated dilation of coronary stents is unknown at present.
- When multiple stents are required, if placement results in stent to stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.
- Polymer coated implants may cause inflammable and/or prothrombotic reactions. The treating physician should balance the medical advantage of a biostable polymer-coated Paclitaxel-Eluting Stent against the potential risk.

## Drug interactions

The amount of Paclitaxel contained in one stent creates its minimum concentration; therefore drug interactions may be neglected.

## Application method

### Precaution Prior to Use

Carefully inspect the sterile package before opening. Do not use APOLLO Coronary Intravascular Stent after the "Expire" date. If the integrity of the sterile package has been compromised prior to the "Expire" date (e.g., damage of the package), contact manufacturer QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the Russian Federation -

SC "Lancet". Contact addresses and phone numbers are given in section "Warranty/Liability" of these instructions for Use. Do not use if any defects are noted in the device or package.

**NOTE:** If at any time during use the proximal shaft for catheter delivery has been bent or kinked, do not continue to use the device.

**Devices and materials required (not included in APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange package)** (All devices and materials are disposable: do not reuse and/or resterilize them):

- Guiding catheter(s) (inner lumen min. 1.42 mm) with the appropriate configuration for the coronary artery to be treated;
- 10 or 20 ml syringe;
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/ 2.39mm;
- Guide wire 0.014 inch/ 0.36 mm;
- 60% contrast medium diluted 1:1;
- Sterile physiological saline solution;
- Inflation device and a three-way stopcock for balloon inflation;
- Introducer.

## Preparation for Use

### Unboxing Removal (Step Action)

1. Carefully remove APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange from its protective tubing for preparation of the implantation system. Do not bend or kink hypotube of the dilatation catheter of implantation system during removal.
2. Remove the transparent cover (cap) from the dilatation catheter from the tip and carefully retrieve the protective guide wire. Holding the dilatation catheter of implantation system with one hand on the proximal part from the stent, with the other hand carefully remove the transparent cover (cap) of the stent distally.

**NOTE:** Care should be taken not to kink or bend the shaft of the dilatation catheter upon application.

### Stent System Guide wire Lumen Flush (Step Action)

1. Flush stent system guide wire lumen with heparinized saline using a syringe.
2. Verify that the stent is positioned between the proximal and distal radiopaque balloon markers. Check for bends, kinks and other damage. Do not use if any defects are noted.

### Preparation of the dilatation catheter of implantation system (Step Action)

Prepare inflation device/syringe with diluted contrast medium.

Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypo tube when connecting to inflation device/syringe.

With tip down, orient stent implantation system vertically.

Open stopcock to stent implantation system; pull negative for 15 seconds; release to neutral for contrast fill.

Close stopcock to stent implantation system; purge inflation device/syringe of all air.

Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use device.

If a syringe was used, attach a prepared inflation device to stopcock.

Open stopcock to the dilatation catheter of stent implantation system.

Leave on neutral.

### **Stent Delivery Procedure (Step Action)**

1. Prepare the vascular access site according to standard percutaneous transluminal coronary angioplasty (PTCA) practice.
2. Predilate the lesion/vessel with appropriate diameter balloon having a ratio of 1:1 with the diameter of the vessel.
3. Maintain neutral pressure on inflation device attached to the dilatation catheter of stent implantation system.
4. Put the dilatation catheter of stent implantation system onto the proximal part of guide wire, while maintaining the distal part of guide wire position across target lesion.
5. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
6. Ensure that the guiding catheter is securely placed in the lesion area, then insert the dilatation catheter of implantation system in the coronary artery. Carefully advance the stent implantation system into the hub of the guiding catheter, keeping the hypo tube of the dilatation catheter of stent implantation system straight.

### **TE**

If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance indicate a problem and use of excessive force may result in stent damage or stent dislodgement from dilatation catheter of stent implantation system. Maintain guide wire placement across the lesion and move the stent implantation system and guiding catheter as a single unit.

7. Advance the stent implantation system over the guide wire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point for correct positioning. The inside edges of the marker bands indicate both the

stent edges and balloon shoulders inflated. If the position of the stent is not optimal, do not deploy the stent; in such case the stent should be carefully repositioned or removed.

3. Sufficiently tighten the rotating hemostatic valve. Stent is now ready to be deployed.

### **Stent Deployment Procedure (Step Action)**

1. Inflate the balloon of catheter of stent implantation system to a minimum pressure of 9 bar, (900 kPa) (deployment pressure). Higher pressures may be necessary to expand the stent to optimize stent apposition against the arterial wall. The pressure of the balloon of catheter of stent implantation system must not exceed rated burst pressure of 18 bar (1800 kPa) for stent diameters 2.00 to 3.50 mm and 16 bar (1600 kPa) for stent diameters at 4.00 mm.
2. Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
3. Deflate balloon by pulling negative on inflation device until balloon is fully deflated.
4. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent.  
Optimal expansion requires that the stent be in full contact with the artery wall. All efforts should be taken to assure that the stent is not underdilated. Stent expansion should be performed with fluoroscopic visualization.
5. If stent sizing/apposition requires optimization, readvance the balloon of the catheter of stent implantation system, or another dilatation catheter of the appropriate size, to the stented area using standard angioplasty techniques.
6. Inflate the balloon to the desired pressure while observing by X-ray imaging. Deflate the balloon (see Balloon Compliance Chart supplied with device).
7. Reconfirm stent position and angiographic result. Ensure that the stent is positioned between the proximal and distal radiopaque markers of the balloon of stent implantation system. Repeat inflations until the desired result is achieved.

### **Stent Implantation System Removal Procedure (Step Action)**

1. Ensure balloon of the catheter is fully deflated.
2. Fully open rotating hemostatic valve.
3. While maintaining guide wire position and negative pressure on inflation device, withdraw implantation system.

### **Precautions when handling APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Before using, check the product "Expire" date.

- Verify visually that the stent is located between the proximal and distal balloon markers of the implantation system.
- Coronary Intravascular Stent Implantation System/Rapid Exchange is designed for use as a unit. The stent is not to be removed from its delivery balloon; it may result in damaging and/or embolization of the stent. The stent is not designed to be crimped onto another balloon.
- Special care must be taken not to handle or in any way disrupt the stent position on the implantation system. This is most important during catheter removal from packaging, placement over guide wire, and advancement through hemostasis valve adaptor and guiding catheter hub.
- Excessive manipulation, e.g. rolling the mounted stent, may result in damaging the coating or dislodgement of the stent from the implantation system.
- Avoid handling the stent on the delivery hands without using protective gloves.
- In the event the stent is not deployed on the implantation system, please contact QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the RF CJSC "Lancet" to clarify the return procedure. Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use.
- Use only the appropriate solution of contrast medium. Do not use air or any gas medium to inflate the balloon.
- Stent contact with any fluid prior to placement is not recommended as there is a possibility of Paclitaxel release. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

## **Precautions when implanting and removing APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Do not pre-inflate balloon prior to stent deployment other than as directed.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- The target lesion must be sufficiently predilated prior to stent implantation.
- Do not expand the stent if it is not properly positioned in the vessel.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the Balloon Compliance Chart fixed to product label. Use of pressures higher than specified on product label may result in a ruptured balloon and potential intimal damage and dissection. The vessel should be pre-dilated with appropriate diameter balloon having a 1:1 ratio with the vessel diameter.

- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.
- If unusual resistance is felt at any time during lesion access before stent implantation, the stent implantation system and the guiding catheter should be removed as a single unit.
- Do not attempt to pull an unexpanded stent back into the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.

**NOTE: The following steps should be executed under direct visualization using X-ray imaging.**

- Maintain guide wire placement across the lesion during the entire removal process. Carefully pull back the stent implantation system until the proximal radiopaque balloon marker of the stent implantation system is aligned with the distal tip of the guiding catheter.
- The stent implantation system and guiding catheter should be pulled back until the tip of the guiding catheter is just distal to the arterial sheath, allowing the guiding catheter to straighten. Carefully retract the stent implantation system into the guiding catheter and remove the stent implantation system and guiding catheter from the patient as a single unit while leaving the guide wire across the target lesion.

Failure to follow these steps and/or applying excessive force to the stent implantation system can potentially result in stent damage, stent dislodgement from the balloon and/or damage to the implantation system.

## **Precautions after implantation of APOLLO Coronary Intravascular Stent**

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, or a dilatation catheter to avoid disrupting the stent geometry or coating.
- Do not perform Magnetic Resonance Imaging (MRI) scan on patient's post-stent implantation until the stent has been completely endothelialized (90 days) to minimize the potential for migration. The stent may cause artefacts in MRI scans due to distortion of the magnetic field.
- Prescribe an antiplatelet therapy for a period of 6 months to reduce the risk of stent thrombosis.

## **Reuse precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

## **Package content:**

1 APOLLO Coronary Intravascular Stent Implantation System/Rapid Exchange in dispenser with balloon compliance scale packed in sealed sterile plastic package. Package with device, instructions for use and desiccant placed in a cardboard box.

## **Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.
2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

## **Disposal instructions**

After use, dispose of product and packaging in accordance with common requirements to disposal of waste of medical devices.

## **Shelf Life**

Shelf life – 2 years from the date of sterilization. Do not use after expiry.

## **Warranty/Liability**

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the Instruction of Use. The warning notices being contained in this Instruction of Use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.

Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Schenstrasse 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, Building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

#### Graphical symbols for labeling primary and secondary packages

| Symbol | Description                                                            |
|--------|------------------------------------------------------------------------|
|        | Sterilized with ethylene oxide gas                                     |
|        | Lot number                                                             |
|        | Manufacturer                                                           |
|        | Each device is suitable for one (1) time use only, reuse is PROHIBITED |
|        | Expires                                                                |
|        | Keep in a cool, dry and dark place                                     |
|        | Read instructions prior to use                                         |
|        | Hydrophilic coating                                                    |

|                                                                                   |                                   |
|-----------------------------------------------------------------------------------|-----------------------------------|
|  | CE Marking                        |
|  | Logo:<br>IK /INTELLIGENT KINETICS |

**Manufacturer**  
**Med Innovative Medizinprodukte GmbH**  
 chstrasse 16  
 23 Winsen  
 many  
 ne No.: +49 4171 6578 0  
 +49 4171 6578 11

## INSTRUCTIONS FOR USE

### ABARIS Coronary Intravascular Stent Implantation System (Rapid Exchange)

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout the instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

ABARIS Coronary Intravascular Stent Implantation System/ Rapid Exchange is a flexible balloon expandable stent mounted on a rapid exchange (RX) stent implantation system. The RX stent implantation system is designed for safe insertion of the stent in the coronary vessels and for stent deployment in the desired area by means of expansion of the balloon on which the stent is mounted.

The RX Stent Implantation System consists of a dilatation catheter with a balloon on the distal end. The stent is mounted on the balloon. The balloon of stent implantation system expands to the nominal diameter at the nominal pressure.

ABARIS Coronary Intravascular Stent Implantation System/ Rapid Exchange has the following features:

- Stent is made from stainless steel (316 LVM), carbon ions impregnated up to 0.5  $\mu\text{m}$  deep, with biodegradable Sirolimus-eluting polymeric coating. Sirolimus content - 2.0  $\mu\text{g}/\text{mm}^2$  ( $\pm 10\%$ ). Sirolimus share – up to 33% (mass). Coating thickness - 5  $\mu\text{m}$ . Polymeric coating - RESOMER RG 504 (INN: PLGA (*Poly lactide 50%, Polyglycolid 50%*)).
- The balloon of the dilatation catheter of stent implantation system is made from nylon Grilamid L25.
- Balloon diameter – from 2 to 4 mm; balloon length from 12 to 40 mm.
- Tip of the dilatation catheter and flexible distal shaft are made from High-Density polyethylene (HDPE) Pebax 7233, and MEVOPUR Black colorant is added to the tip of the dilatation catheter.
- The dilatation catheter of stent implantation system is fitted with a three-layer inner tube made from High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant.
- In order to ensure best sliding and stent implantation to the target area, the distal shaft of the dilatation catheter of stent implantation system is fitted with two-layer hydrophilic coating. Base coat - DSM Primer 41002, Top coat - DSM Comfort Coat Topcoat 43003.
- Hypo tube of the dilatation catheter of stent implantation system consists of two parts: 1) (PTFE) Polytetrafluoroethylene (Xylan) coated stainless steel 304 L and 2) uncoated stainless steel 304 L.

- Kinking protection cover of the dilatation catheter of stent implantation system is made from Thermoplastic Urethane Texin RxT85A with Blue tioxide TR-27 colorant.
- The cannula Luer Lock is made from Polycarbonate (PC) Lexan HPS 7. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint.
- The balloon of the dilatation catheter of stent implantation system has two radiopaque markers (Platinum – 90%, Iridium– 10%), used for precise positioning of the balloon.
- The balloon of stent implantation system can be used for post-dilatation purposes.
- The dilatation catheter of stent implantation system is fitted with two shaft markers made from stainless steel 304 L, designed to determine the position of the dilatation catheter implantation system in relation to the guiding catheter.
- Balloon nominal pressure – 9 bar.
- Burst pressure – 18 bar (for balloon diameter from 2 to 3.50 mm); 16 bar (for balloon diameter 4.00 mm).
- Dilatation catheter of stent implantation system is compatible with guide wire diameter 0.014”.
- The tip of dilatation catheter of stent implantation system is put on the protective guide wire made from stainless steel WKN 1.4301 (tip of protective guide wire from Polytetrafluoroethylene (PTFE) Teflon with Blue basic colorant), length 360, diameter 0.38 mm, passing on the inside of the distal part of the catheter and comes outside through the port.
- In order to protect the stent and the balloon from damage, on the side of the tip the dilatation catheter is fitted with a transparent cover (cap) made from Polytetrafluoroethylene (PTFE) Teflon, which should be removed before the stenting procedure along with retrievable protective guide wire.

**Table 1. Dimensions of ABARIS Coronary Intravascular Stent, dimensions of the balloon of dilatation catheter of stent implantation system**

| Stent length, mm | Balloon length, mm | Usable balloon length, ± 0.5 mm | Inner distance of the markers, ± 0.5 mm | Strut thickness, mm | Strut breadth, mm |
|------------------|--------------------|---------------------------------|-----------------------------------------|---------------------|-------------------|
| 10               | 12                 | 12,50                           | 11,30                                   | 0,110 ± 0,010       | 0,120 ± 0,010     |
| 14               | 15                 | 15,50                           | 14,30                                   |                     |                   |
| 18               | 20                 | 20,50                           | 19,30                                   |                     |                   |
| 24               | 25                 | 25,50                           | 24,30                                   |                     |                   |
| 28               | 30                 | 30,50                           | 29,30                                   |                     |                   |
| 34               | 35                 | 35,50                           | 34,30                                   |                     |                   |
| 38               | 40                 | 41,50                           | 40,30                                   |                     |                   |
| 5                | 10                 | 12,50                           | 11,30                                   | 0,110 ± 0,010       | 0,120 ± 0,010     |
|                  | 14                 | 15,50                           | 14,30                                   |                     |                   |

|   |    |    |       |       |               |                |
|---|----|----|-------|-------|---------------|----------------|
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,110 ± 0,010 | 0,120 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 5 | 10 | 12 | 12,50 | 11,30 | 0,110 ± 0,010 | 0,120 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,110 ± 0,010 | 0,120 ± 0,010  |
|   | 14 | 15 | 15,50 | 14,30 |               |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 5 | 10 | 12 | 12,50 | 11,30 | 0,120 ± 0,010 | 0, 130 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 | 0,115 ± 0,010 |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,120 ± 0,010 | 0, 130 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 | 0,115 ± 0,010 |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |
| 0 | 10 | 12 | 12,50 | 11,30 | 0,120 ± 0,010 | 0, 130 ± 0,010 |
|   | 14 | 16 | 16,50 | 15,30 | 0,115 ± 0,010 |                |
|   | 18 | 20 | 20,50 | 19,30 |               |                |
|   | 24 | 25 | 25,50 | 24,30 |               |                |
|   | 28 | 30 | 30,50 | 29,30 |               |                |
|   | 34 | 35 | 35,50 | 34,30 |               |                |
|   | 38 | 40 | 41,50 | 40,30 |               |                |

## Table 2. Technical characteristics of dilatation catheter of stent implantation system

| Parameter                                                                                  | Value                                             |
|--------------------------------------------------------------------------------------------|---------------------------------------------------|
| Available length of dilatation catheter of stent implantation system, mm                   | 1415 ± 15                                         |
| Length of dilatation catheter of stent implantation system, mm                             | 3,5 ± 0,5                                         |
| Distal part of catheter, mm                                                                | 237±15                                            |
| Flexible distal shaft (Transition zone from hypo tube to outer tubing of the catheter), mm | 105±15                                            |
| Hypo tube, mm                                                                              | 1070±15                                           |
| Proximal shaft diameter of the dilatation catheter of stent implantation system, F / mm    | Max. 2,1 F (0,70 mm)                              |
| Distal shaft diameter of the dilatation catheter of stent implantation system, F / mm      | Max. 2,6 F (0,86 mm)                              |
| Proximal Shaft marker Distance of dilatation catheter, mm                                  | 985 ± 15                                          |
| Distal Shaft marker Distance of dilatation catheter, mm                                    | 880 ± 15                                          |
| Working protection cover, mm                                                               | Diameter – 1,10<br>Length - 35<br>Thickness – 0,2 |
| Annula Luer Lock, mm                                                                       | Length - 41<br>Inner diameter from 0,7 to 1,25    |
| Distal balloon length of the dilatation catheter of stent implantation system, mm          | Depends on stent size*                            |
| Available balloon length of the dilatation catheter of stent implantation system, mm       | Depends on stent size*                            |
| Marker distance of the markers                                                             | Depends on stent size*                            |
| tapered Tip Profile of the dilatation catheter of stent implantation system, mm            | 0,45                                              |
| Proximal entrance profile of the dilatation catheter of stent implantation system, mm      | 0,54                                              |

### Indications for use

ABARIS Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for treatment of coronary arteries when treating coronary vessel constriction due to advancing arterial sclerotic disease, not requiring open surgery.

ABARIS Coronary Intravascular Stent Implantation System/ Rapid Exchange is designed for improving the lumen of native coronary arteries (de novo) or vein grafts:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;

- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

Drug eluting of ABARIS coronary intravascular stent significantly helps reduce:

- restenosis frequency after percutaneous coronary intervention,
- revascularization frequency of target lesion,
- frequency of late loss of vessel lumen after stenting procedure.

Implantation of ABARIS Drug-Eluting Coronary Intravascular Stent is also indicated for treatment of acute or threatening closure in patients with failed interventional therapy.

The treated lesion length should be less than the nominal stent length (stent length - from 10 to 38 mm) with reference vessel diameters from 2.0 to 4.0 mm.

## Contraindications

Implantation of ABARIS intravascular coronary stent implantation system/rapid exchange is contraindicated for the patients having the following conditions:

- Hypersensitivity to Sirolimus.
- Allergic reaction to stainless steel.
- Allergic reaction to PLGA (Polylactide 50%, Polyglycolid 50%).
- Severe reaction to contrast agents.
- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.
- In-stent restenosis.
- Myocardial infarction < 72 hours.
- Stenting of saphenous vein grafts.
- Unprotected left main coronary artery.
- Total occlusion of target vessel.
- Heavily calcified lesions.
- Lesions involving arterial segments with highly tortuous anatomy.
- Lesions involving a bifurcation. Stent implantation in bifurcation area.
- Left ventricular ejection fraction < 30 %.
- Cardiogenic shock.
- Presence of definite or probable intraluminal thrombus.
- Any patients judged to have a lesion which may interfere with proper stent deployment.

- Pregnancy and breast-feeding.

## **essential Adverse Events**

- Abrupt stent closure.
- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.
- Stenocardia.
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).
- Arteriovenous fistula.
- Cardiac tamponade.
- Cardiogenic shock.
- Death.
- Emboli, distal (air, tissue or thrombotic emboli).
- Emergent coronary artery bypass surgery.
- Heart failure.
- Hematoma.
- Haemorrhage.
- Hypotension/Hypertension.
- Infectious complications.
- Infectious inflammation.
- Myocardial ischemia.
- Coronary artery dissection, perforation or rupture.
- Pericardial effusion.
- Pseudo aneurysm, femoral.
- Renal failure.
- Restenosis of stented segment.
- Shock/Pulmonary edema.
- Coronary artery spasm.
- Stent embolization.
- Stent thrombosis/occlusion.
- Stroke/cerebrovascular accident/transient ischemic attack.
- Total occlusion of coronary artery or shunt.
- Vessel trauma requiring surgical repair or reintervention.
- Histologic changes in vessel wall, including inflammation, cellular damage or necrosis.

Allergic reaction to stainless steel.

Allergic reaction to Sirolimus, to PLGA (Polylactide 50%, Polyglycolid 50%).

Alopecia.

Myalgia/Arthralgia.

Gastrointestinal symptoms.

Peripheral neuropathy.

Hematologic dyscrasia (anaemia, leukopenia, neutropenia, thrombocytopenia).

Hepatic enzyme changes.

## arnings

Application of ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange carries an associated risk of subacute thrombosis, vascular complications such as bleeding, hematoma or pseudoaneurysm. Antiplatelet therapy must be prescribed for all patients for a period of 6 months post procedure.

ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange must not be implanted in patients with allergic reactions to stainless steel, PLGA (Polylactide 50%, Polyglycolid 50%) and/or hypersensitivity to Sirolimus.

Implantation of ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange should be performed only by physicians who have received appropriate training.

Implantation of ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.

Direct stenting without prior dilation of the lesion has not been studied using this product.

Implantation of ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange may lead to development of subsequent restenosis, which may require reintervention (particularly, insertion and expansion of a dilatation catheter in arterial segment containing the stent). The long-term outcome following repeated dilation of coronary stents is unknown at present.

When multiple stents are required, if placement results in stent to stent contact, stent materials should be of similar composition to avoid the possibility of dissimilar metal corrosion.

## g interactions

amount of Sirolimus contained in one stent creates its minimum concentration; therefore drug actions may be neglected.

## Application method

### Precaution Prior to Use

Carefully inspect the sterile package before opening. Do not use ABARIS Coronary Intravascular Stent Implantation System after the "Expire" date. If the integrity of the sterile package has been compromised prior to the "Expire" date (e.g., damage of the package), contact manufacturer QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the Russian Federation - CJSC "Medica". Contact addresses and phone numbers are given in section "Warranty/Liability" of these instructions for Use. Do not use if any defects are noted in the device or package.

**NOTE:** If at any time during use the proximal shaft for catheter delivery has been bent or kinked, do not continue to use the device.

**Devices and materials required (not included in ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange package)** (All devices and materials are disposable; do not reuse or resterilize them):

- Guiding catheter(s) (inner lumen min. 1.42 mm) with the appropriate configuration for the coronary artery to be treated;
- 10 or 20 ml syringe;
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/ 2.39mm;
- Guide wire 0.014 inch/ 0.36 mm;
- 60% contrast medium diluted 1:1;
- Sterile physiological saline solution;
- Inflation device and a three-way stopcock for balloon inflation;
- Introducer.

### Preparation for Use

#### Stent Packaging Removal (Step Action)

1. Carefully remove ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange from its protective tubing for preparation of the implantation system. Do not bend or kink hypo tube of the dilatation catheter of implantation system during removal.
2. Remove the transparent cover (cap) from the dilatation catheter from the tip and carefully retrieve the protective guide wire. Holding the dilatation catheter of implantation system with one hand on the proximal part from the stent, with the other hand carefully remove the transparent cover (cap) of the stent distally.

**NOTE:** Care should be taken not to kink or bend the shaft of dilatation catheter upon application.

## **Stent System Guide wire Lumen Flush (Step Action)**

1. Flush stent system guide wire lumen with heparinized saline using a syringe.
2. Verify that the stent is positioned between the proximal and distal radiopaque balloon markers.  
Check for bends, kinks and other damage. Do not use if any defects are noted.

## **Preparation of dilatation catheter of implantation system (Step Action)**

- Prepare inflation device/syringe with diluted contrast medium.
- Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypo tube when connecting to inflation device/syringe.
- With tip down, orient stent implantation system vertically.
- Open stopcock to stent implantation system; pull negative for 15 seconds; release to neutral for contrast fill.
- Close stopcock to stent implantation system; purge inflation device/syringe of all air.
- Repeat steps 4 through 6 until all air is expelled. If bubbles persist, do not use device.
- If a syringe was used, attach a prepared inflation device to stopcock.
- Open stopcock to dilatation catheter of stent implantation system.
- Leave on neutral.

## **Stent Delivery Procedure (Step Action)**

1. Prepare the vascular access site according to standard percutaneous transluminal coronary angioplasty (PTCA) practice.
2. Predilate the lesion/vessel with appropriate diameter balloon having a ratio of 1:1 with the diameter of the vessel.
3. Maintain neutral pressure on inflation device attached to dilatation catheter of stent implantation system.
4. Put the dilatation catheter of stent implantation system onto the proximal part of guide wire, while maintaining the distal part of guide wire position across target lesion.
5. Fully open rotating hemostatic valve to allow for easy passage of the stent and prevent damage to the stent.
6. Ensure that the guiding catheter is securely placed in the lesion area, then insert the dilatation catheter of implantation system in the coronary artery. Carefully advance the stent implantation system into the hub of the guiding catheter, keeping the hypo tube of the dilatation catheter of stent implantation system straight.

If unusual resistance is felt before the stent exits the guiding catheter, do not force passage. Resistance indicate a problem and use of excessive force may result in stent damage or stent dislodgement from dilatation catheter of stent implantation system. Maintain guide wire placement across the lesion and move the stent implantation system and guiding catheter as a single unit.

7. Advance the stent implantation system over the guide wire to target lesion under direct fluoroscopic visualization. Utilize the proximal and distal radiopaque balloon markers as a reference point for correct positioning. The inside edges of the marker bands indicate both the stent edges and balloon shoulders inflated. If the position of the stent is not optimal, do not deploy the stent; in such case the stent should be carefully repositioned or removed.
8. Sufficiently tighten the rotating hemostatic valve. Stent is now ready to be deployed.

### Stent Deployment Procedure (Step Action)

1. Inflate the balloon of catheter of stent implantation system to a minimum pressure of 9 bar, (900 kPa) (deployment pressure). Higher pressures may be necessary to expand the stent to optimize stent apposition against the arterial wall. The pressure of the balloon of catheter of stent implantation system must not exceed rated burst pressure of 18 bar (1800 kPa) for stent diameters 2.00 to 3.50 mm and 16 bar (1600 kPa) for stent diameters at 4.00 mm.
2. Maintain inflation pressure for 15-30 seconds for full expansion of the stent.
3. Deflate balloon by pulling negative on inflation device until balloon is fully deflated.
4. Confirm stent position and deployment using standard angiographic techniques. For optimal results, the entire stenosed arterial segment should be covered by the stent.  
Optimal expansion requires that the stent be in full contact with the artery wall. All efforts should be taken to assure that the stent is not underdilated. Stent expansion should be performed with fluoroscopic visualization.
5. If stent sizing/apposition requires optimization, readvance the balloon of the catheter of stent implantation system, or another dilatation catheter of the appropriate size, to the stented area using standard angioplasty techniques.
6. Inflate the balloon to the desired pressure while observing by X-ray imaging. Deflate the balloon (see Balloon Compliance Chart supplied with device).
7. Reconfirm stent position and angiographic result. Ensure that the stent is positioned between the proximal and distal radiopaque markers of the balloon of stent implantation system. Repeat inflations until the desired result is achieved.

## **Implantation System Removal Procedure (Step Action)**

- Ensure balloon of the catheter is fully deflated.
- Fully open rotating hemostatic valve.
- While maintaining guide wire position and negative pressure on inflation device, withdraw implantation system.

## **Precautions when handling ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Before using, check the product "Expire" date.
- Verify visually that the stent is located between the proximal and distal balloon markers of the implantation system.
- Coronary Intravascular Stent Implantation System/Rapid Exchange is designed for use as a unit. The stent is not to be removed from its delivery balloon, it may result in damaging and/or embolization of the stent. The stent is not designed to be crimped onto another balloon.
- Special care must be taken not to handle or in any way disrupt the stent position on the implantation system. This is most important during catheter removal from packaging, placement over guide wire, and advancement through hemostasis valve adaptor and guiding catheter hub.
- Excessive manipulation, e.g. rolling the mounted stent, may result in damaging the coating or dislodgement of the stent from the implantation system.
- Avoid handling the stent on the delivery hands without using protective gloves.
- In the event the stent is not deployed on the implantation system, please contact QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the RF CJSC "Lancet" to clarify the return procedure. Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use.
- Use only the appropriate solution of contrast medium. Do not use air or any gas medium to inflate the balloon.
- Stent contact with any fluid prior to placement is not recommended as there is a possibility of Sirolimus release. However, if it is absolutely necessary to flush the stent with sterile/isotonic saline, contact time should be limited (1 minute maximum).

## **Precautions when implanting and removing ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange**

- Do not pre-inflate balloon prior to stent deployment other than as directed.
- Implanting a stent may lead to dissection of the vessel distal and/or proximal to the stented portion, and may cause acute closure of the vessel requiring additional intervention (e.g., CABG, further dilation, placement of additional stents, or other).
- The target lesion must be sufficiently predilated prior to stent implantation.
- Do not expand the stent if it is not properly positioned in the vessel.
- Balloon pressures should be monitored during inflation. Do not exceed rated burst pressure as indicated on the Balloon Compliance Chart fixed to product label. Use of pressures higher than specified on product label may result in a ruptured balloon and potential intimal damage and dissection. The vessel should be pre-dilated with appropriate diameter balloon having a 1:1 ratio with the vessel diameter.
- Stent retrieval methods (use of additional wires, snares and/or forceps) may result in additional trauma to the vascular site. Complications can include bleeding, hematoma or pseudoaneurysm.
- If unusual resistance is felt at any time during lesion access before stent implantation, the stent implantation system and the guiding catheter should be removed as a single unit.
- Do not attempt to pull an unexpanded stent back into the guiding catheter while engaged in the coronary arteries, as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.
- An unexpanded stent should be introduced into the coronary arteries one time only. An unexpanded stent should not be subsequently moved in and out through the distal end of the guiding catheter as stent damage or stent dislodgement from the balloon of catheter of stent implantation system may occur.

### **NOTE: The following steps should be executed under direct visualization using X-ray imaging.**

- Maintain guide wire placement across the lesion during the entire removal process. Carefully pull back the stent implantation system until the proximal radiopaque balloon marker of the stent implantation system is aligned with the distal tip of the guiding catheter.
- The stent implantation system and guiding catheter should be pulled back until the tip of the guiding catheter is just distal to the arterial sheath, allowing the guiding catheter to straighten. Carefully retract the stent implantation system into the guiding catheter and remove the stent implantation system and guiding catheter from the patient as a single unit while leaving the guide wire across the target lesion.

Failure to follow these steps and/or applying excessive force to the stent implantation system can potentially result in stent damage, stent dislodgement from the balloon and/or damage to the implantation system.

**Precautions after implantation of ABARIS Coronary Intravascular Stent**

- Care must be exercised when crossing a newly deployed stent with an intravascular ultrasound (IVUS) catheter, a coronary guide wire, or a dilatation catheter to avoid disrupting the stent geometry or coating.
- Do not perform Magnetic Resonance Imaging (MRI) scan on patient's post-stent implantation until the stent has been completely endothelialized (90 days) to minimize the potential for migration. The stent may cause artefacts in MRI scans due to distortion of the magnetic field.
- Prescribe an antiplatelet therapy for a period of 6 months to reduce the risk of stent thrombosis.

**Use precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Package content:** 1 ABARIS Coronary Intravascular Stent Implantation System/Rapid Exchange in dispenser with balloon compliance scale packed in sealed sterile plastic package. Package with device, instructions for use and desiccant placed in a cardboard box.

**Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.
2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

**Disposal instructions**

After use, dispose of product and packaging in accordance with common requirements to disposal of waste of medical devices.

## If Life

Shelf life – 1 year from the date of sterilization. Do not use after expiry.


## arranty/Liability

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the instructions of use. The warning notices being contained in this instruction of use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.

Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Schwanenstraße 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, Building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

## Graphical symbols for labeling primary and secondary packages

| Symbol                                                                              | Description                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Sterilized with ethylene oxide gas                                     |
|  | Lot number                                                             |
|  | Manufacturer                                                           |
|  | Each device is suitable for one (1) time use only, reuse is PROHIBITED |
|  | Do not resterilize                                                     |
|  | Expires                                                                |

|                                                                                     |                                    |
|-------------------------------------------------------------------------------------|------------------------------------|
|     | Keep in a cool, dry and dark place |
|    | Read instructions prior to use     |
|    | Hydrophilic coating                |
|    | CE Marking                         |
|  | Logo:<br>IK /INTELLIGENT KINETICS  |

**Manufacturer**  
**AliMed Innovative Medizinprodukte GmbH**  
 Schwanenbachstrasse 16  
 30634 Winsen  
 Germany  
 Phone No.: +49 4171 6578 0  
 Fax No.: +49 4171 6578 11

## INSTRUCTIONS FOR USE

### PEGASUS PTCA Balloon Catheter

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

PEGASUS PTCA Balloon Catheter is a coronary dilatation catheter designed to enlarge the lumen of coronary arteries or venous shunts (with the diameter from 1.5 to 5.0 mm) narrowed by atherosclerotic lesions.

PEGASUS PTCA Balloon Catheter is a double lumen catheter with hydrophilic coating and a balloon at the distal tip of the catheter. One lumen of the catheter is adapted to deliver contrast agent to the balloon for balloon inflation. The other lumen of distal shaft allows using the guide wire to facilitate the advancement of balloon catheter to and through the stenosis.

The desired balloon size depends on the atherosclerotic lesion size. The balloon should cover more than 50% of lesion.

- The balloon of the PEGASUS Catheter is made from nylon Grilamid L25.
- Balloon diameter – from 1.5 to 5 mm; balloon length from 10 to 40 mm.
- Tip of the PEGASUS Balloon Catheter and flexible distal shaft are made from High-Density polyethylene (HDPE) Pebax 7233, and MEVOPUR Black colorant is added to the tip of the PEGASUS Balloon Catheter.
- The PEGASUS Balloon Catheter is fitted with a three-layer inner tube made from High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant.
- In order to ensure best sliding and balloon delivery to the target area, the distal shaft of the PEGASUS Balloon Catheter is fitted with two-layer hydrophilic coating. Base coat - DSM Primer 41002, Top coat - DSM Comfort Coat Topcoat 43003.
- Hypo tube of the PEGASUS Balloon Catheter consists of two parts: 1) (PTFE) Polytetrafluoroethylene (Xylan) coated stainless steel 304 L and 2) uncoated stainless steel 304 L.
- Kinking protection cover of the PEGASUS Balloon Catheter is made from Thermoplastic Urethane Texin RxT85A with Blue tioxide TR-27 colorant.
- Cannula Luer Lock is made from Polycarbonate (PC) Lexan HPS 7. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint.
- The balloon of the PEGASUS Catheter has two radiopaque markers (Platinum – 90%, Iridium – 10%), used for precise positioning of the balloon.

- The PEGASUS Balloon Catheter is fitted with two shaft markers made from stainless steel 304 L, designed to determine the position of the balloon catheter in relation to the guiding catheter tip.
- Balloon nominal pressure – 9 bar.
- Burst pressure – 18 bar (for balloon diameter from 1.50 to 3.50 mm); 16 bar (for balloon diameter from 4.00 to 4.50 mm); 14 bar (for balloon diameter 5.00 mm).
- The PEGASUS Balloon Catheter is compatible with guide wire diameter 0.014”.
- The tip of the PEGASUS Balloon Catheter is put on the protective guide wire made from stainless steel WKN 1.4301, length 325 mm, diameter 0.38 mm, passing on the inside of the distal part of the catheter and comes outside through the port.
- In order to protect the balloon from damage, PEGASUS catheter is equipped with an orange cover (cap) made of Polytetrafluoroethylene (PTFE) Teflon with Orange classic colorant.

**Table 1. Dimensions of PEGASUS Balloon Catheter**

| Balloon diameter, mm | Balloon length, mm | Usable balloon length, mm<br>± 0.5 mm | Inner distance of the markers, mm<br>± 0.5 mm |
|----------------------|--------------------|---------------------------------------|-----------------------------------------------|
| 1,5                  | 10                 | 9,50                                  | 8,3                                           |
|                      | 11                 | 11,50                                 | 10,3                                          |
|                      | 12                 | 12,50                                 | 11,3                                          |
|                      | 15                 | 15,50                                 | 14,3                                          |
|                      | 16                 | 16,50                                 | 15,3                                          |
|                      | 20                 | 20,50                                 | 19,3                                          |
| 2                    | 10                 | 9,50                                  | 8,3                                           |
|                      | 11                 | 11,50                                 | 10,3                                          |
|                      | 12                 | 12,50                                 | 11,3                                          |
|                      | 15                 | 15,50                                 | 14,3                                          |
|                      | 16                 | 16,50                                 | 15,3                                          |
|                      | 20                 | 20,50                                 | 19,3                                          |
|                      | 25                 | 25,50                                 | 24,3                                          |
|                      | 30                 | 30,50                                 | 29,3                                          |
|                      | 35                 | 35,50                                 | 34,3                                          |
| 2,25                 | 10                 | 9,50                                  | 8,3                                           |
|                      | 11                 | 11,50                                 | 10,3                                          |
|                      | 12                 | 12,50                                 | 11,3                                          |
|                      | 15                 | 15,50                                 | 14,3                                          |
|                      | 16                 | 16,50                                 | 15,3                                          |
|                      | 20                 | 20,50                                 | 19,3                                          |
|                      | 25                 | 25,50                                 | 24,3                                          |
|                      | 30                 | 30,50                                 | 29,3                                          |
|                      | 35                 | 35,50                                 | 34,3                                          |
| 2,50                 | 10                 | 9,50                                  | 8,3                                           |
|                      | 11                 | 11,50                                 | 10,3                                          |
|                      | 12                 | 12,50                                 | 11,3                                          |

|      |    |       |      |
|------|----|-------|------|
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 2,75 | 10 | 9,50  | 8,3  |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 3,00 | 10 | 9,50  | 8,3  |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 3,25 | 10 | 9,50  | 8,3  |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 3,50 | 10 | 9,50  | 8,3  |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 4,00 | 10 | 9,50  | 8,3  |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |

|      |    |       |      |
|------|----|-------|------|
|      | 15 | 15,50 | 14,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |
| 4,50 | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
| 5,00 | 40 | 41,50 | 40,3 |
|      | 11 | 11,50 | 10,3 |
|      | 12 | 12,50 | 11,3 |
|      | 16 | 16,50 | 15,3 |
|      | 20 | 20,50 | 19,3 |
|      | 25 | 25,50 | 24,3 |
|      | 30 | 30,50 | 29,3 |
|      | 35 | 35,50 | 34,3 |
|      | 40 | 41,50 | 40,3 |

**Table 2. Technical characteristics of PEGASUS Balloon Catheter**

| Parameter                                                                                              | Value                                             |
|--------------------------------------------------------------------------------------------------------|---------------------------------------------------|
| Usable length of PEGASUS Balloon Catheter, mm                                                          | 1415 ± 15                                         |
| Tip of PEGASUS Balloon Catheter, mm                                                                    | 3,5 ± 0,5                                         |
| Distal part of PEGASUS Balloon Catheter, mm                                                            | 237±15                                            |
| Flexible distal shaft (Transition zone from hypo tube to inner tubing of PEGASUS Balloon Catheter), mm | 105±15                                            |
| Hypo tube, mm                                                                                          | 1070±15                                           |
| Proximal shaft diameter of PEGASUS Balloon Catheter, F / mm                                            | Max. 2,1 F (0,70 mm)                              |
| Distal shaft diameter of PEGASUS Balloon Catheter, F / mm                                              | Max. 2,6 F (0,86 mm)                              |
| Tip - Proximal Shaft marker Distance of PEGASUS Balloon Catheter, mm                                   | 985 ± 15                                          |
| Tip - Distal Shaft marker Distance of PEGASUS Balloon Catheter, mm                                     | 880 ± 15                                          |
| Kinking protection cover, mm                                                                           | Diameter – 1,10<br>Length - 35<br>Thickness – 0,2 |
| Cannula Luer Lock, mm                                                                                  | Length – 41<br>Inner diameter from 0,7 to 1,25    |
| Total balloon length of PEGASUS Balloon Catheter                                                       | Depends on stent size*                            |
| Usable balloon length of PEGASUS Balloon Catheter                                                      | Depends on stent size*                            |

|                                                      |                        |
|------------------------------------------------------|------------------------|
| Inner distance of the markers                        | Depends on stent size* |
| Tapered Tip Profile of PEGASUS Balloon Catheter, mm  | 0,45                   |
| Tip entrance profile of PEGASUS Balloon Catheter, mm | 0,54                   |

### Indications for use

PEGASUS Balloon Catheter is designed to enlarge the lumen of coronary vessels with the diameter from 1.5 to 5.0 mm in the treatment of progressive atherosclerotic stenosis that does not require open surgery.

PEGASUS Balloon Catheter is designed for direct deployment during dilatation and for pre- and post-dilatation of coronary arteries.

PEGASUS Balloon Catheter is designed to enlarge the lumen of coronary vessels, in particular before stenting procedures:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;
- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

The treated lesion length should be less than the nominal balloon length (balloon length – from 10 to 40 mm) with reference vessel diameters from 1.5 to 5.0 mm.

### Contraindications

PEGASUS Balloon Catheter is contraindicated for the patients having the following conditions:

- Severe reaction to contrast agents.
- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.
- Unprotected left main coronary artery.
- Pregnancy and breast-feeding.
- Patients with a history of gastrointestinal hemorrhage.
- Surgical treatment two weeks before implantation, act of delivery, blood vessel biopsy and puncture of vessel not subject to compression.
- Patients with recent acute cerebrovascular event (ACE).
- Diabetic hemorrhagic retinopathy and other conditions associated with coagulation defects.
- Coronary artery spasm with no significant stenosis.

## Possible adverse events

- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.
- Stenocardia.
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).
- Arteriovenous fistula.
- Death.
- Vessel dissection.
- Hematoma.
- Hemorrhage.
- Hypotension/Hypertension.
- Infection complications.
- Infectious inflammation.
- Myocardial ischemia.
- Coronary artery dissection, perforation or rupture.
- Restenosis of stented segment.
- Coronary artery spasm.
- Total occlusion of coronary artery or shunt.
- Vessel trauma requiring surgical repair or reintervention.
- Coronary embolism.

## Warnings

- In order to decrease the risk of vessel injury, the diameter of inflated balloon should approximate to the diameter of the vessel just proximal and distal to the stenosis.
- All procedures with PEGASUS Balloon Catheter shall be done under direct X-ray visualization.
- PEGASUS Balloon Catheter shall be withdrawn only after the balloon is deflated completely by applying negative pressure.
- The pressure of the balloon of PEGASUS Balloon Catheter must not exceed rated burst pressure.
- Before PEGASUS Balloon Catheter is advanced into the vessel, evacuate all air from the balloon.
- Use only the appropriate solution of contrast medium. Do not use air or any gas medium to inflate the balloon.
- If at any time during use the proximal shaft for PEGASUS Balloon Catheter has been bent or kinked, do not continue to use the device.

- PEGASUS Balloon Catheter may be applied only by certified medical professionals with proper training and skills in angiography and percutaneous transluminal coronary angioplasty (PTCA) of coronary vessels. During intervention, a cardiac surgical team should be on standby.
- PTCA should only be performed at hospitals where emergency coronary artery bypass graft surgery can be readily performed.
- In the patients not eligible for coronary artery bypass surgery, PTCA requires all factors to be considered thoroughly, in particular, hemodynamic support during PTCA procedure.
- In the patients with diabetes mellitus, the risk of possible complications is higher.

## **Application method**

### **Inspection Prior to Use**

Thoroughly inspect the sterile package before opening. Do not use PEGASUS Balloon Catheter after the "Expire" date. Check PEGASUS Balloon Catheter for bends, torsions or other defects. If the integrity of the sterile package has been compromised prior to the product "Expire" date (e.g., damage of the package), contact manufacturer QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the Russian Federation - CJSC "Lancet". Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use. Do not use if any defects are noted in the device or package.

Devices and materials required for angioplasty and stenting procedures (not included in PEGASUS Balloon Catheter package) (All devices and materials are disposable: do not reuse and/or resterilize):

- Guiding catheter(s) (inner lumen min. 1.42 mm) with the appropriate configuration for the coronary artery to be treated;
- 10 or 20 ml syringe;
- Rotating hemostatic valve(s) with a minimum diameter of 0.096 inch/ 2.39mm;
- Guide wire 0.014 inch/ 0.36 mm;
- 60% contrast medium diluted 1:1;
- Sterile physiological saline solution;
- Inflation device and a three-way stopcock for balloon inflation;
- Introducer.

### **Preparation for Use**

#### **Preparation and Removal (Step Action)**

1. Carefully remove PEGASUS Balloon Catheter from its protective tubing for preparation of the delivery system. Do not bend or kink hypo tube of PEGASUS Balloon Catheter during removal.
2. Remove the transparent cover (cap) from the balloon catheter from the tip and carefully retrieve the protective guide wire.

**NOTE: Care should be taken not to kink or bend the shaft of balloon catheter upon application.**

### **Flushing of guide wire lumen of the catheter (Step-by-step instruction)**

1. Flush the guide wire lumen of the catheter by injecting heparinized saline from syringe connected to a port. The heparinized saline shall flow from the guide wire lumen.
2. Be sure that radiopaque markers in the proximal and distal parts of PEGASUS catheter balloon are positioned in the target lesion site. Check for bends, kinks and other damage. Do not use if any defects are noted.

### **Preparation of PEGASUS Balloon Catheter (Step Action)**

Prepare inflation device/syringe with diluted contrast medium.

Attach inflation device/syringe to stopcock; attach to inflation port. Do not bend the hypo tube when connecting to inflation device/syringe.

With tip down, orient PEGASUS Balloon Catheter vertically.

Open stopcock to PEGASUS Balloon Catheter; pull negative for 15 seconds; release to neutral for contrast fill.

Disconnect the syringe or balloon inflation device from the port of PEGASUS catheter.

Close stopcock to PEGASUS Balloon Catheter; purge inflation device/syringe of all air.

Reconnect the syringe or balloon inflation device to the balloon inflation port of PEGASUS catheter.

Maintain negative pressure in the balloon until air stops entering the catheter.

Release pressure to zero.

Disconnect the syringe and connect the balloon inflation device to the balloon inflation port of PEGASUS catheter, without introducing air into the system.

### **Application of PEGASUS Balloon Catheter for coronary artery dilatation**

Advance the guide wire through the hemostatic valve of introducer per standard procedure. Then advance guide wire carefully through the guiding catheter. Then withdraw the guide wire introducer, if used.

Attach a torque device to the guide wire.

Under X-ray visualization, advance the guide wire to the target vessel, then across the stenosis.

When backloading the balloon catheter onto the guide wire on the distal end, be sure that the guide wire exits the catheter. Advance PEGASUS balloon catheter through the guide wire up to the hemostatic valve. Loosen the knurled knob of hemostatic valve. Advance the catheter through the hemostatic valve introducer, while maintaining the guide wire position, and tighten the knurled knob of hemostatic valve.

To facilitate catheter insertion, the balloon should be fully deflated to negative pressure.

NOTE: When backloading the catheter onto the guide wire, the catheter shall be supported by hands. In advancing the catheter into the guide wire, support it with one hand and grasp firmly the proximal part of catheter.

Tighten the knurled knob of hemostatic valve around the PEGASUS Balloon Catheter to form a seal without inhibiting the movement of catheter. This will allow continuous recording of proximal coronary artery pressure.

Advance the balloon catheter until the appropriate proximal radiopaque marker aligns the knurled knob of hemostatic valve. Such aligning of balloon marker to the knurled knob of hemostatic valve indicates that the distal tip of balloon catheter has reached the distal tip of guiding catheter.

Advance the balloon catheter over the guide wire to the stenosis. Inflate the balloon at a low pressure (1-2.0 MPa), to confirm that the balloon is positioned correctly.

Deflate the balloon and perform PTCA per standard procedure. Between inflations, maintain negative pressure on the balloon.

Withdraw the deflated balloon catheter and guide wire from the guiding catheter through the hemostatic valve. Tighten the knurled knob of hemostatic valve.

## Precautions

1. Before using, check the product "Expire" date.
2. Examine PEGASUS Balloon Catheter to ensure that it is not damaged, its dimensions are suitable for the specific procedure for which it is to be used.
3. During the procedure, appropriate anticoagulant, anti-platelet and coronary vasodilator therapy must be provided to the patients. Such therapy should be continued for a period of time determined by a physician.

4. Do not use PEGASUS Balloon Catheter without gloves.
5. When opening and closing the hemostatic valve of introducer, and when withdrawing the catheter, the differences in catheter shaft diameter shall be taken into account.
6. The hemostatic valve should not be closed so tight that it restricts the flow of contrast agent into and out of the balloon of catheter. Provided that the hemostatic valve should be closed tightly enough to prevent the blood leakage around the catheter.
7. When using the dual wire technique, a dual hemostatic valve should be used. Care shall be taken to avoid entanglement of wires during insertion, torque or withdrawal. When using the dual wire technique, guide wires should not be rotated more than 180° in either direction.
8. It is recommended that one wire is completely withdrawn from the patient blood system before removing additional equipment.
9. This catheter is not designed for distal artery pressure recording.

#### **Reuse precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing and/or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Package content:** 1 PEGASUS Balloon Catheter in dispenser with balloon compliance scale packed in sealed sterile plastic package. Package with device, instructions for use and desiccant placed in a cardboard box.

#### **Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.
2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

#### **Disposal instructions**

After use, dispose of product and packaging in accordance with common requirements to disposal of waste of medical devices.

## Shelf Life

Shelf life - 2 years from the date of sterilization. Do not use after expiry.



## Warranty/Liability

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the Instruction of Use. The warning notices being contained in this Instruction of Use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.



Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Schwanenstrasse 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, Building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

## Graphical symbols for labeling primary and secondary packages

| Symbol                                                                              | Description                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|  | Sterilized with ethylene oxide gas                                     |
|  | Lot number                                                             |
|  | Manufacturer                                                           |
|  | Each device is suitable for one (1) time use only, reuse is PROHIBITED |
|  | Expires                                                                |

|                                                                                     |                                    |
|-------------------------------------------------------------------------------------|------------------------------------|
|    | Keep in a cool, dry and dark place |
|    | Read instructions prior to use     |
|    | Hydrophilic coating                |
|    | CE Marking                         |
|  | Logo:<br>IK /INTELLIGENT KINETICS  |

**Manufacturer**  
**aliMed Innovative Medizinprodukte GmbH**  
 Schstrasse 16  
 23 Winsen  
 rmany  
 one No.: +49 4171 6578 0  
 : +49 4171 6578 11 T

## INSTRUCTIONS FOR USE

### POLLUX Aspiration Catheter (Set)

Carefully read all instructions prior to use. Observe all warnings and precautions noted throughout these instructions; otherwise, complications may occur. For single patient use only. Use the product before expiration of the „Best-Before-Date” mentioned on the packing.

**STERILE:** Sterilized with ethylene oxide. Do not use if the sterile package is open or damaged.

#### Description

POLLUX Aspiration Catheter (Set) is designed to remove thrombotic burden by percutaneous aspiration from coronary arteries before stenting in patients suffering thrombosis.

POLLUX Aspiration Catheter (Set) is composed of:

1. POLLUX Aspiration Catheter compatible with 5F, 6F or 7F Guiding Catheter.
2. 2 Syringes, 30 ml each.
3. Cell strainer.
4. Stopcock with extension tubing.

**POLLUX Aspiration Catheter** is a dual lumen rapid exchange (RX) catheter compatible with 5F, 6F or 7F guide catheters. One lumen (inner tube) of the catheter is designed to aspirate thrombotic burden or soft emboli from coronary arteries. A single lumen tube of the catheter is made of High-Density polyethylene (HDPE) Pebax 7233. The second lumen of the catheter (a three-layer inner tube housed inside the rapid exchange segment) is designed to advance a guide wire to the place of aspiration. The three-layer inner tube of rapid exchange segment is made of High-Density polyethylene (HDPE) Pebax 7233 / Low-Density polyethylene (LDPE) Plexar PX 3080 / High-Density polyethylene (HDPE) Petrothene-LM 6007 with Reflex blue C colorant. Tip of POLLUX Aspiration Catheter is made of High-Density polyethylene (HDPE) Pebax 7233 with Reflex blue C colorant.

The Luer Lock cannula of POLLUX Aspiration Catheter is made of Polycarbonate (PC) Makrolon 2458. Signs are made on Luer Lock cannula wings with No-Tox Black vinyl paint. The Luer Lock cannula is supported by a tube for proximal stiffness from Polyester Vention with Blue (Opaque) colorant.

The stainless steel 304 V stylet loaded in POLLUX Aspiration Catheter assists the advancement of catheter. Luer Lock cover is made of Polycarbonate (PC) Makrolon 2458 with Orange C59 colorant. The stylet allows distributing the effort required for advancing catheter. The stylet shall be removed from POLLUX Aspiration Catheter before aspiration.

The tip of POLLUX Aspiration Catheter is put on the protective guide wire made of WKN 1.4301 stainless steel with the length of 325 mm and diameter of 0.38 mm, passing on the inside of the distal part of the catheter and coming outside from the port. The protective guide wire is designed to protect the distal part of POLLUX Aspiration Catheter from breaks. The protective guide wire shall be removed from the catheter before aspiration.

In order to protect the proximal part of catheter from break/kink, the catheter is equipped with the protective tip made of thermoplastic urethane Texin RxT85A with White tioxide TR-23 colorant.

**Table 1. Technical characteristics of POLLUX Aspiration Catheter**

| Parameter                                                                        |          | Value               |                     |                     |
|----------------------------------------------------------------------------------|----------|---------------------|---------------------|---------------------|
| Diameter of compatible guide catheter, F / ° / mm                                |          | 5 F / 0.058° / 1,47 | 6 F / 0.070 °/ 1,78 | 7 F / 0.081 °/ 2,06 |
| Usable length of POLLUX Aspiration Catheter, mm                                  |          | 1410 ± 15           | 1365 ± 15           | 1410 ± 15           |
| Diameter of guide wire, °/ mm                                                    |          | 0,014°/0.36         | 0,014°/0.36         | 0,014°/0.36         |
| Length of rapid exchange segment, mm                                             |          | 75 ± 3              | 30 ± 3              | 75 ± 3              |
| Distance from the radiopaque marker to the tip of POLLUX Aspiration Catheter, mm |          | 3,5 ± 1,1           | 3,5 ± 1,1           | 3,5 ± 1,1           |
| Length of aspiration opening, mm                                                 |          | 4,8 ± 1,1           | 4,8 ± 1,1           | 4,8 ± 1,1           |
| Shaft diameter, mm                                                               | Distal   | 1,067 ± 0,025       | 1,37 ± 0,05         | 1,55 ± 0,03         |
|                                                                                  | Central  |                     |                     |                     |
|                                                                                  | Proximal | 1,12 ± 0,025        | 1,42 ± 0,05         |                     |
| Length of the tube securing the proximal stiffness, mm                           |          | 300±25              | 300±25              | 300±25              |
| Length of protective tip against break/kink, mm                                  |          | 38 ± 0,05           | 38 ± 0,05           | 38 ± 0,05           |
| Minimal diameter of Luer Lock cannula, mm                                        |          | 0,7                 | 0,7                 | 0,7                 |
| Outer diameter of Luer Lock outlet fitting, mm                                   |          | 7,82                | 7,82                | 7,82                |
| Inner diameter of Luer Lock outlet fitting, mm                                   |          | 4,3                 | 4,3                 | 4,3                 |
| Length of the thread on the Luer Lock thread fitting, mm                         |          | 6,6                 | 6,6                 | 6,6                 |
| Length of Luer Lock outlet fitting, mm                                           |          | 7,87                | 7,87                | 7,87                |
| Total length of Luer Lock outlet and thread fittings, mm                         |          | 15,74               | 15,74               | 15,74               |

|                                                                             |         |         |         |
|-----------------------------------------------------------------------------|---------|---------|---------|
| Total length of Luer Lock fitting, mm                                       | 17,78   | 17,78   | 17,78   |
| Outer diameter of Luer Lock output fitting, mm                              | 3,9     | 3,9     | 3,9     |
| Inner diameter of Luer Lock thread fitting, mm                              | 7,9     | 7,9     | 7,9     |
| Distal shaft, mm                                                            | 360 ±15 | 360 ±15 | 360 ±15 |
| Distance from the port to the edge of hydrophilic coating, mm               | 10      | 10      | 10      |
| Distance from the end of hydrophilic coating to the end of distal shaft, mm | 10      | 10      | 10      |

POLLUX Aspiration Catheter has a smooth atraumatic tip with 1 radiopaque marker (Platimun - 90%, Iridium - 10%) on the distal end for catheter precise positioning. To ensure smooth and efficient delivery, the distal shaft of POLLUX catheter has a two-layer hydrophilic coating, the basecoat is made of DSM Primer 41002, and the topcoat is made of DSM Comfort Coat Topcoat 43003.

**Syringe** has a polycarbonate barrel (PC Lexan HPS 7) with a 6-position locking device (plunger). Signs are made on barrel with No-Tox Black vinyl paint. The plunger tip is made of silicone Dow Corning QP1 with Black 108 colorant, the plunger is made of ABS (Acrylonitrile butadiene styrene) Lustran ABS 248 White. The syringe is designed for the aspiration of thrombotic burden. If the aspirated volumes exceed the syringe volume, the first syringe shall be disconnected after the stopcock on extension tubing is closed and another syringe shall be connected. For this purpose, the Set of POLLUX Aspiration Catheter includes two syringes of 30 ml each.

**Table 2. Technical characteristics of syringe, 30 ml**

| Parameter                                 | Value, cm |
|-------------------------------------------|-----------|
| Syringe length from tip to plunger flange | 12,42     |
| Length of syringe barrel, including tip   | 9,9       |
| Width of plunger stopper                  | 2,54      |
| Length of syringe handle                  | 5,71      |
| Width of syringe barrel                   | 2,64      |
| Thickness of syringe handle               | 0,18      |

**Cell strainer** is designed for isolating blood clots from the aspirated contents. The cell strainer is a polypropylene frame (PP Purell RP271 with Blue CC 324 C colorant) with the attached 40 µm nylon mesh (Nylon BD Falcon).

**Table 3. Technical characteristics of cell strainer**

| Parameter      | Value, mm |
|----------------|-----------|
| Height         | 14,60     |
| Lower diameter | 25,2      |
| Upper diameter | 25,7      |
| Handle length  | 11,1      |

**Stopcock with extension tubing** is designed to connect the syringe to the catheter. The tubing is made of PVC Xtra-Flex HS-914, stopcock is made of Polyoxymethylene / Polyacetal (POM) Delrin with Green H35 colorant, adaptor for connection with syringe and adaptor for connection with Luer Lock cannula - Polycarbonate (PC) Makrolon 2458 are made of Polycarbonate (PC) Makrolon 2458.

**Table 4. Technical characteristics of Stopcock with extension tubing**

| Parameter                  | Value, cm  |
|----------------------------|------------|
| Length, stopcock inclusive | 19 ±1      |
| Outer diameter             | 0,125±0,03 |
| Inner diameter             | 0,06±0,02  |

**Indications for use**

POLLUX Aspiration Catheter (Set) is designed for the removal of thrombotic burden or soft emboli from coronary arteries, in particular before stenting procedures:

- in patients with symptomatic ischemic disease due to stenotic lesions of native coronary vessels;
- in patients with symptomatic ischemic disease due to occlusive lesions of native coronary vessels;
- in patients with symptomatic ischemic disease after coronary artery bypass grafting when treating stenotic lesions of coronary artery bypasses;
- for treatment of acute or threatening occlusion of coronary arteries resulted from failed percutaneous coronary intervention.

**Contraindications**

POLLUX Aspiration Catheter (Set) is contraindicated for the patients having the following conditions:

- Patients in whom anti-platelet and/or anticoagulant therapy is contraindicated.

- Unprotected left main coronary artery.
- Pregnancy and breast-feeding.
- Cerebral thrombosis.
- Venous thrombosis.
- Surgical treatment two weeks before implantation, act of delivery, blood vessel biopsy and puncture of vessel not subject to compression.
- Lesions involving a bifurcation.
- Heavily calcified lesions.

### **Potential Adverse Events**

- Acute myocardial infarction.
- Allergic reaction to anticoagulant and/or anti-platelet therapy and/or contrast medium and/or anesthetic agents.
- Stenocardia.
- Arrhythmias, including ventricular fibrillation (VF) and ventricular tachycardia (VT).
- Arteriovenous fistula.
- Cardiac tamponade.
- Death.
- Emboli, distal (air, tissue or thrombotic emboli).
- Emergent coronary artery bypass surgery.
- Hematoma.
- Haemorrhage.
- Hypotension/Hypertension.
- Infectious complications.
- Infectious inflammation.
- Myocardial ischemia.
- Coronary artery dissection, perforation or rupture.
- Pseudo aneurysm, femoral.
- Renal failure.
- Restenosis of coronary artery.
- Shock/Pulmonary edema.
- Coronary artery spasm.
- Stroke/cerebro vascular accident/transient ischemic attack.
- Total occlusion of coronary artery / Thrombosis.

- Vessel trauma requiring surgical repair.
- Histologic changes in vessel wall, including inflammation, cellular damage or necrosis.

## **Warnings and Precautions**

1. Before using, check the product “Expire” date.
2. POLLUX Aspiration Catheter (Set) may be applied only by certified medical professionals with proper training and skills in angiography and percutaneous transluminal coronary angioplasty (PTCA) of coronary vessels. During intervention, a cardiac surgical team should be on standby.
3. Do not expose to organic solvents or oil-based contrast agents. The use of such contrast media may lead to damage of the aspiration catheter.
4. Check for bends, kinks and other damage. Do not use if any defects are noted.
5. Prior to use, ensure that the diameter and length of the aspiration catheter are correctly matched to the guide catheter and guide wire.
6. Check that all connections of the set are secure so that air is not introduced into the syringe or extension tubing through the catheter during aspiration.
7. When aspirating thrombus, adjust the aspiration volume considering the diameter of vessel from which the thrombus is sucked and the blood flow into the target lesion. (If these two parameters are too small compared to the aspiration volume of this catheter, negative pressure may arise at the distal tip of catheter.)
8. If the guide catheter is wedged in the entrance of the coronary artery or if there is stenosis upstream of the target lesion which is subject to aspiration procedure, the procedure shall be performed cautiously as the vessel or implanted stent may be deformed by a negative pressure arising during aspiration. [If the blood flow at the target lesion where the aspiration of thrombus is performed is too small compared to the aspiration volume of this catheter, negative pressure may arise at the distal tip of catheter.]
9. When using the catheter for aspirating thrombus from the carotid artery, a distal tip protection device is mandatory to occlude temporarily the distal side of the target lesion.
10. When advancing this catheter into vessels, the guide catheter shall be infused with heparinized saline to prevent coagulation.
11. The catheter shall be prepared for use by removing residual air from the catheter and syringe.
12. The catheter can only be inserted with the use of a guide wire, otherwise the damage or proliferation of the vascular wall may occur.

13. The guide wire shall be completely advanced to the end of vessel containing the lesion to be aspirated, otherwise it may loosen from the guide wire lumen and cause damage to vascular wall or proliferation of vessel.
14. If the guide wire is displaced from the guide wire lumen during procedure, remove the catheter and reinsert the guide wire.
15. During aspiration procedure, the position of catheter and thrombus shall be verified by X-ray visualization. [During the insertion of this catheter into the area of thrombus, the catheter may move the thrombus and cause it to disperse. In addition, attachment of the distal end of the catheter onto the vascular wall during aspiration may result in the damage to the vascular wall.]
16. If strong resistance is experienced during the procedure, the cause of such resistance shall be found and measures shall be taken before proceeding. Otherwise the vessels may be damaged or catheter shaft may break and remain inside the patient.
17. If this catheter is advanced to the area where the stent is placed, insertion or pulling out shall be performed slowly to avoid the distal end of the catheter becoming stuck on a stent strut. The interface between the catheter and the guide wire may catch the stent strut, which could potentially damage the stent and the catheter.  
If any resistance is encountered, do not advance or remove the catheter with effort. Check whether the catheter is caught and if so try to remove it from the stent strut by adjusting the position of catheter and the guide wire.
18. Do not advance the catheter to the distal part of site where the venous filter is implanted: the strut of venous filter may damage the catheter or break the catheter shaft leaving debris inside the body.
19. If broken or kinked, the shaft shall be replaced. [If such defect is ignored and catheter continues to be used, aspiration may fail or the catheter shaft may be damaged and remain in the patient].
20. If a greater resistance is encountered during insertion, advancement or pulling out of this catheter, it shall be verified that the guide wire is not tangled. If so (taking into consideration the short guide wire lumen of this catheter), the guide wire may wind around the catheter shaft. In addition, while drawing this catheter back into the guide catheter inside vessels, the wide-angle separation between the catheter shaft and guide may occur. Under this circumstance, a forced withdrawal of the catheter may lead to damage to the catheter or guide wire.
21. Due to the lack of conductivity of twist forces, the catheter shaft may be twisted. If twisted, the catheter shaft may be damaged and then remain in the patient.

22. The aspiration volume shall be decided by the physician individually for each lesion and adjusted via the aspiration syringe.
23. Air and blood sucked into syringe from thrombus shall not be injected from the catheter into the body.
24. This catheter shall not be inserted or pulled swiftly out from lesions with highly tortuous vessels, bifurcated lesions or calcified lesions. [In this case, the shaft at the distal end may kink or break resulting in the damage to the vessel.]
25. Precautions should be taken to prevent any damage to the catheter by other equipment (such as scalpels, blades or scissors). Do not use a damaged catheter.
26. When using this catheter, the temperature, blood pressure, pulse, and respiration of patients shall be monitored. In case of any abnormality, the procedure shall be stopped or appropriate measures taken based on physician's judgment.
27. Do not use the catheter for infusing diagnostic, embolic or therapeutic materials into blood vessels, since the catheter has not been designed for such use.
28. Do not use the syringe, extension tubing or stopcock inside the human body.
29. Refer to the instructions supplied with other interventional devices to be used in conjunction with this set for their intended use, contraindications, precautions and instructions for use.
30. If blood does not fill syringe during aspiration and catheter blockage is suspected, do not flush the catheter while in patient.
31. Do not perform high pressure contrast injections into the guide catheter while using with the aspiration catheter inside. This may damage the aspiration catheter and make it difficult to remove it from the guide catheter.
32. If flow of aspirated contents is restricted, do not attempt to flush the catheter lumen while the catheter is still inside the patient. It may result in serious complications, such as intravascular thrombus delivery, thromboembolic event and/or death of patient. The catheter shall be removed and flushed outside the patient or a new catheter shall be used.
33. During using of catheter, appropriate drug therapy (anticoagulants, vasolidators, etc.) shall be administered to the patient according to standard protocols.

## **Application method**

### **Inspection Prior to Use**

Carefully inspect the sterile package before opening. Do not use POLLUX Aspiration Catheter (Set) after the "Expire" date. If the integrity of the sterile package has been compromised prior to the product "Expire" date (e.g., damage of the package), contact

manufacturer QualiMed Innovative Medizinprodukte GmbH, Germany or its authorized representative in the Russian Federation - CJSC "Lancet". Contact addresses and phone numbers are given in section "Warranty/Liability" of these Instructions for Use. Do not use if any defects are noted in the device or package.

**NOTE: If at any time during use the proximal shaft for POLLUX Aspiration Catheter has been bent or kinked, do not continue to use the device.**

**Devices and materials required, necessary for aspiration (not supplied in the package of POLLUX Aspiration Catheter (Set))** (All devices and materials are disposable: do not reuse and/or resterilize them):

- Y-connector.
- Guide wire.
- Guide catheter.
- Heparinized saline.
- Balloon device for distal protection (optional).
- Flushing needle.
- Introducer.


## Preparation for Use

### Packaging Removal (Step Action)

1. Aseptically remove this catheter set from the hermetically sealed sterile plastic package. Remove POLLUX Aspiration Catheter from dispenser.
2. Remove protection wire from the distal part of the aspiration catheter.

### Flushing of Aspiration and Guide Wire Lumens of POLLUX Aspiration Catheter (Step-by-Step Instruction)

- a) Draw 10 ml of heparinized saline into the syringe supplied in set and expel any air from it.
- b) Attach the syringe tip onto the adaptor of stopcock extension tubing and attach the adaptor of extension tubing to the Luer Lock connector of POLLUX Aspiration Catheter (as shown on the picture below)

- 1 – Cannula Luer Lock of POLLUX Aspiration Catheter
- 2 – Luer Lock connector of POLLUX Aspiration Catheter
- 3 - Stopcock of extension tubing
- 4 – Syringe

- c) Air inside the aspiration lumen of POLLUX Aspiration Catheter and in the syringe shall be replaced with the full volume of heparinized saline inside the syringe. In the case of incomplete air replacement, procedures 3a) to 3c) shall be repeated, and after flushing the stopcock of extension tubing shall be closed.
- d) Attach the flushing needle to another prepared syringe filled with heparinized saline. Fill the guide wire lumen at the distal tip of catheter with saline.

**Application of POLLUX Aspiration Catheter (Set) for Aspiration**

**Procedures for Insertion of Catheter into the Body**

- 1) Prior to use, check if POLLUX Aspiration Catheter is compatible to the guide catheter and the guide wire according to the following table:

| Compatible guide catheter diameter, F"/mm | Guide wire diameter,"/mm |
|-------------------------------------------|--------------------------|
| 5 F / 0.058" / 1.47 mm                    | 0,014"/0.36 mm           |
| 6 F / 0.070" / 1.78 mm                    | 0,014"/0.36 mm           |
| 7 F /0.081"/2.06 mm                       | 0,014"/0.36 mm           |

- 2) Insert the guide catheter and after the installation of Y-connector advance it more distally from the vascular lesion. If distal protection balloon is used, advance the balloon to the end of the lesion.
- 3) Insert the proximal end of the guide wire into the distal part of the corresponding catheter lumen.
- 4) Insert POLLUX Aspiration Catheter through Y-connector under the X-ray visualization to make the radiopaque marker at the distal tip of the catheter reach the target site.  
To prevent shaft twisting, hold POLLUX catheter 2- 3 cm from Y-connector, insert slowly the catheter into the introducer entrance until the catheter reaches the target site.
- 5) Moderately lock the Y-connector to prevent blood loss and resistance to catheter manipulations during the aspiration procedure.
- 6) Disconnect the adaptor of extension tubing (connected to the syringe) from LuerLock connector of stylet, remove stylet from POLLUX catheter.

7) Reattach the adaptor of extension tubing to Luer Lock cannula of POLLUX Aspiration Catheter.

8) After verifying the location of the catheter under the X-ray visualization, pull the plunger of the syringe to the desired volume. After locking the plunger, reduce the pressure inside the syringe.

9) Open the stopcock of extension tubing and start aspiration:

**If no protection balloon is used:**

Open the stopcock of extension tubing and advance slowly POLLUX catheter to the thrombus for aspiration. Then draw back POLLUX catheter to the proximal end of lesion.

**If the balloon device is used for distal protection and prevention of embolization with thrombotic material,** in accordance with instruction for use of protection balloon, temporally occlude the target lesion. Then open the stopcock and advance slowly POLLUX catheter through the whole area of thrombus for aspiration.

Then close the stopcock of extension tubing and remove the syringe. If urgent recanalization is required, deflate immediately the distal protection balloon according to its instruction for use.

10) If another aspiration is required, reattach the syringe to the stopcock of extension tubing after removing the aspirated material from the syringe and repeat procedures 8 and 9.

**Removal of POLLUX Aspiration Catheter (Set)**

- 1) Ensure the aspiration of the thrombus is completed;
- 2) Ensure the stopcock of the extension tubing is closed;
- 3) Remove slowly POLLUX Aspiration Catheter;
- 4) Ensure no thrombus remains inside the Y-connector. If necessary, remove thrombus material from Y-connector.
- 5) If the reinsertion of POLLUX Aspiration Catheter is necessary, the insertion should be performed again for the priming procedure after the insertion of the stylet into POLLUX catheter.

**Reuse precaution statement**

For single patient use only. Do not reuse, reprocess and/or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.

**Package content:** 1 POLLUX Aspiration Catheter in dispenser, 2 syringes, 1 mesh, 1 extension tubing with stopcock attached to the backing card and packed hermetically into the sterile plastic bag. Package with devices, instructions for use and desiccant placed in a card board box.

### **Recommended transportation, storage and usage conditions**

1. Transport at temperature 20 - 25°C, relative humidity 70%.
2. Keep in a cool, dark and dry place, at temperature 20 - 25°C, relative humidity 70%.
3. Usage conditions- at temperature 20 - 25°C, relative humidity 70%.
4. See section "Graphical symbols for labeling primary and secondary packages".
5. Use immediately after the sterile package is opened.

### **Disposal instructions**

After use, dispose of product and packaging in accordance with common requirements to disposal of waste of medical devices.

### **Shelf Life**

Shelf life – 3 years from the date of sterilization. Do not use after expiry.



### **Warranty/Liability**

Carefully read instructions of use before proceeding. QualiMed Innovative Medizinprodukte GmbH warranty exists only under the provision that the use of the product has been carried out in accordance with the Instruction of Use. The warning notices being contained in this Instruction of Use have to be considered particularly as an essential part of warranty granted by QualiMed Innovative Medizinprodukte GmbH.

Use the product before expiration of the „Best-Before-Date” mentioned on the packing. QualiMed Innovative Medizinprodukte GmbH grants a warranty on the products up to the mentioned “Best-Before-Date”.

If the integrity of the sterile package has been compromised, if any other defects in the package or the device have been identified, please contact QualiMed Innovative Medizinprodukte GmbH, Boschstrasse 16, 21423 Winsen, Germany, phone No.: +49 4171 6578 0, fax: +49 4171 6578 11 or its authorized representative Closed Joint-Stock Company "Lancet" (CJSC "Lancet"), 17, Otkrytoe Shosse, building 1, Moscow, Russia, 107143, phone No.: + 7 (495) 646-56-65, fax: + 7 (495) 646-56-66.

Graphical symbols for labeling primary and secondary packages

| Symbol                                                                              | Description                                                            |
|-------------------------------------------------------------------------------------|------------------------------------------------------------------------|
|    | Sterilized with ethylene oxide gas                                     |
|    | Lot number                                                             |
|    | Manufacturer                                                           |
|    | Each device is suitable for one (1) time use only, reuse is PROHIBITED |
|    | Do not resterilize                                                     |
|   | Expires                                                                |
|  | Do not use if package is compromised                                   |
|  | Keep in a cool, dry and dark place                                     |
|  | Read instructions prior to use                                         |
|  | Hydrophilic coating                                                    |

|                                                                                   |                                   |
|-----------------------------------------------------------------------------------|-----------------------------------|
|                                                                                   |                                   |
|  | CE Marking                        |
|  | Logo:<br>IK /INTELLIGENT KINETICS |

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
Unseren Chief Risk Officer (Risiko-Manager):

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N: DE50 240 603 00 420 392 580 0

BIC-/SWIFT-Code: GENODEF1NBU

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I, Michael Lammers, Notary Public, Rathausstr. 60, 21423 Winsen (Luhe), Germany, herewith certify that the overleaf signature attested by me have been served before me by:

Mr. Manfred Gülcher, date of birth: may 04 1961,  
Boschstr. 16, 21423 Winsen (Luhe)  
- personally known to me -

Mr. Gülcher acted by authority. You find the original document enclosed.

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