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
China Council for the Promotion of International Trade is China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of HEBEI XINLE SCI & TECH CO., LTD. on the annexed INSTRUCTION FOR USE is genuine.

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朱广宇
Zhu Guangyu

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(Date: Mar. 16, 2022)

证明书查询网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

<u>INSTRUCTION FOR USE</u>	
Approved Name: BAI QIANG Position: EXPORT MANAGER	Signature, Stamp  河北鑫乐医疗器械科技有限公司 HEBEI XINLE SCITECH CO., LTD

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Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey),
 Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey)
 (see Appendix)

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Name of the tube¹	Cap (plug) color
MINI tube with sodium fluoride and potassium oxalate MINI tube with sodium fluoride and EDTA K2	grey

Potassium oxalate anticoagulant and glucose stabilizer sodium fluoride are contained. Anticoagulant prevents blood from clotting by binding calcium ions, and sodium fluoride inhibits glycolysis by blocking the enzyme enolase. The concentration of glucose in the sample is stable for up to 6 hours at room temperature. The study of glucose concentration should be carried out as soon as possible, since sodium fluoride, inhibiting one of the last stages of glucose destruction, is less effective against the initial stages.

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Attention! Glucose test tubes should be filled to the volume indicated on them, since an excess of anticoagulant can cause hemolysis. The contents of the tube after blood collection should be mixed with extreme caution.

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Purpose: Test tubes are used to determine the concentration of glucose, lactate and glycosylated hemoglobin in plasma or oxalate blood, as well as in tests for glucose tolerance.

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Through the capillary into the mini-test tube with a rubber cap and a cross shaped valve (the minitube should be held horizontally). Wait until the capillary is filled with blood and turn the tube vertically so that the blood can freely enter it from the capillary.

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Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
NAF + Potassium oxalate or with EDTA K2	24 hours	8-10 times	1500 g	10

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Name	Volume	Dimensions	Additive	Cap (plug) color	Capillary	Cap
MINI tube with sodium fluoride and potassium oxalate (GREY)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	sodium fluoride and potassium oxalate	Grey	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
MINI tube with sodium fluoride and EDTA K2 (GREY)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	sodium fluoride and K2EDTA	Grey	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover

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E-mail:lab@int-med.ru

APPENDIX

1. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
2. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
3. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
4. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
5. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
6. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
7. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
8. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
9. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
10. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
11. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
12. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
13. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

14. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
15. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
16. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and potassium oxalate (Grey), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
17. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
18. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
19. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
20. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
21. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
22. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
23. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
24. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
25. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
26. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
27. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary

- blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
28. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 29. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 30. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 31. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
 32. Mini tube non-vacuum for capillary blood (sterile, with sodium fluoride and EDTA K2(Grey), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

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

证明书查询网站 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

<u>INSTRUCTION FOR USE</u>	
Approved Name: BAI QIANG Position: EXPORT MANAGER	Signature, Stamp 

INSTRUCTION FOR USE

Tube non-vacuum, single use, for human capillary blood with additives, with capillaries and without capillaries w in variations for laboratory tests

Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green) Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green)

Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green)

(see Appendix)

RC № RZN XXX from XX.XX.XXXX

DESCRIPTION

Tube non-vacuum, single use, for human capillary blood with additives, with capillaries and without capillaries in variations for laboratory tests – these are cylindrical containers of different heights and diameters, with half-round bottom in volumes from 0.2 to 1 ml, made of PP.

The tube is made of a material that allows you to clearly see the contents of the container when visually inspected. A variety of the products allows you to pick up the suitable tube for specific laboratory research or analysis.

The test tube consists of a body with a cap (plug) into which a capillary can be inserted (for products with a rubber cap with a cross-shaped valve). There are two types of caps (plugs) - rubber with a cross-shaped valve and screw-down. The additive is applied inside the inner wall of the tube.

The capillary is a long thin transparent tube made of polypropylene with a diameter of 2.8 mm. Capillaries are sterile.

The capillaries' length is 53mm.

Type of analyzed sample - Capillary blood sample.

All tubes and capillary are sterile, equipped with a convenient label with a field for recording, which helps to easily identify the tube. Each tube has a calibration mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive)

Name of the tube¹	Cap (plug) color
MINI tube with sodium heparin (NaHe)	green
MINI tube with lithium heparin (LiHe)	
MINI tube with sodium heparin NaHe and separating Gel	

Heparin is an anticoagulant that activates antithrombins and prevents blood clotting, while not distorting the parameters studied. The blood components in a sample with a lithium salt of heparin are stable for up to 6 hours.

Caution: Do not use these tubes for monitoring drugs and determining lithium content.

¹ Further in the text, it is possible to use abbreviated names of versions without specifying volumes, type of lid and packaging and/or other parameters.

IFU non vac (LH/NH/NH+ gel): 01 - 2020

DESIGNATION AND INDICATIONS (purpose)

The product is intended for use in the capture and preservation and / or transportation of human capillary blood for various types of clinical laboratory studies of samples taken.

The product is intended for use in medical institutions, including blood transfusion stations, research institutions. The product is an aid in in vitro diagnostics.

Purpose: Test tubes are used for biochemical and immunological studies of plasma, for studies of heparinized blood and blood cell counts.

CONTRAINDICATIONS AND SIDE EFFECTS

There are no contraindications and side effects at this time.

Precautions. Warnings. CAUTIONS

Precautions:

1. Handle all biological samples and sharp objects (scarificators, lancets) in accordance with the rules adopted at your medical facility.
2. Do not disregard cases of contact with biological specimens (for example, in case of an accidental injection by a lancet), since it can transmit pathogens of dangerous infections (HIV, hepatitis, .).
3. Discard used sharps into a special container for disposal.
4. All liquid reagents are clear. Do not use tubes if the additives have changed color or there is a foreign substance in them.
5. Do not use tubes after the expiration date.

Warnings:

1. Please check the directions for use and designation.
2. Check the tube for damage, contamination, leakage.
3. Disposable medical device. Do not reuse. Dispose after use.
4. Verify that the tube specification and test parameter match.
5. After collecting the blood sample, move the tube carefully. Stir the contents of the tube several times to mix the blood sample with the filler (additive).
6. Avoid cases when collected blood returns to the patient's bloodstream.
7. Be careful when working, avoid contamination.

Cautions and storage conditions (labeled)

- To store at a temperature from +4 ° C to 25 ° C
- Do not use after expiration date.

INSTRUCTION FOR USE

Before the procedure, make sure that the following items are directly accessible:

1. All tubes for capillary blood collection of the required volumes and with the necessary fillers (additives).
2. Capillaries (using tubes with a cross-shaped valve).
3. Lancets or scarifiers for puncture.
4. Personal protective equipment against contact with blood (disposable gloves, overalls, mask).
5. Alcohol wipes for cleaning the puncture site.
6. Dry sterile gauze or dry cotton swabs.
7. Container for recycling or its equivalent.

Technique of capillary blood taking:

1. Disinfect hands and wear disposable gloves.
2. Prepare a mini test tube:

Open the cap of the screw cap tube designed for gravity collection of blood by turning it counterclockwise.

OR

Remove the capillary from the tube and insert into the tube of the second type through the cross-shaped valve on the cap (plug) if there is no capillary, open the cap to collect blood by gravity.

3. Disinfect the puncture site. Prepare the automatic lancet or scarifier according to the manufacturer's instructions.

Attention! Do not use cotton swabs for disinfection, as their fibers can get into the sample and lead to incorrect analysis results (for example, clog the cell count sensor and hemoglobin chamber).

4. Pinch the patient's finger between the thumb, middle and index fingers. OR

For newborns: gently, but firmly take the child by the heel.

5. Make a puncture, remove the first protruding drop.

6. Gently press down on the surrounding tissue so that a drop of blood appears. Without touching the puncture site, collect the dripping blood:

By gravity into a test tube with a cap that turns with the help of a tray ledge. OR

Through the capillary into the mini-test tube with a rubber cap and a cruciform valve (the minitube should be held horizontally). Wait until the capillary is filled with blood and turn the tube vertically so that the blood can freely enter it from the capillary.

Attention! Watch out for the filling level of the minitube. Compliance with the required level of filling allows you to achieve the correct ratio of blood and filler (supplements).

7. After filling the tube:

Close the screw cap of the mini test tube tightly. OR

Carefully remove the capillary from the minitube with a rubber cap and a crossshaped valve (the valve of the stopper closes automatically). The used capillary should be immediately disposed of in a special container or its equivalent.

Attention! If you collect blood in several mini tubes, a new capillary should be used for each.

8. Gently invert the mini-tube 5-7 times to ensure thorough mixing of the filler (additive) and blood. Tap lightly on the bottom of the microfuge to ensure that the filler (additive) and blood are well mixed.

Attention! Do not shake the tubes, as this may cause churning or hemolysis. In tubes with anticoagulants, improper mixing can lead to gluing of platelets, clot formation and incorrect results of the analysis.

9. When all necessary samples have been collected, press a dry sterile swab into the puncture site and wait until the bleeding stops. (If blood was taken from the heel, lift the baby's leg up.)

10. After the blood is drawn, the traces of blood may remain in the recess of the mini-vial cap with the rubber cap and the cross-shaped valve. Take precautions to avoid contact with it.

Centrifuging

Make sure that the mini test tube caps do not rest on the walls of the centrifuge cap, otherwise the cover may come off. See the table below for sample preparation conditions for analysis. The recommended temperature is 15-25 ° C, higher temperatures can have a negative effect on the physical properties of the gel. In centrifuges with horizontal rotors, a more stable gel barrier is formed than in centrifuges with angular rotors.

Preparation of samples for analysis

Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
LH/NH	3 hours	8-10 times	2000 - 3000 g	15
NH+Gel	3 hours before centrifugation and 48 hours after centrifugation	8-10 times	2200 g	15

SYMBOLS WITH DESCRIPTION

	Batch code		Keep dry
	Used by date		Keep away from sunlight
	Manufacturer		Fragile
	Date of manufacturer		Up
	Do not reuse		Layer constraint
	Do not use if pack is damaged		CE marking
	Sterilized by radiation		In vitro diagnostic medical device
	Temperature limit 4-25 °C		

GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.
 Values of the parameters, characteristics (features) of the medical device are presented in the table below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	Volume	Dimensions	Additive	Cap (plug) color	Capillary	Cap type
MINI tube with sodium heparin (NH) (GREEN)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	NH	Green	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
MINI tube with lithium heparin (LH) (GREEN)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	LH	Green	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
MINI tube with sodium heparin NH and separating Gel (GREEN)	0,2 ml, 0,25 ml, 0,5 ml, 0,75 ml	10x45 mm	NH + gel	Green	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover

IFU non vac (LH/NH/NH+ gel): 01 - 2020

MAIN RAW MATERIALS

Tube and capillary PP

Rubber cap and cross-shaped valve Butyl rubber

Screw-down cover PE

STERILIZATION

The products are sterile (tubes and capillaries).

The individually wrapped test tube is sterilely packed together with the capillary.

The tube in a rack and shrink-wrapped is sterile.

The capillaries in the bag, which is attached to the outside of the rack, are also sterile.

Sterilization method - radiation sterilization.

PACKAGE/SCOPE OF DELIVERY

1. Single Use non-vacuum mini tubes for capillary blood with a rubber cap and a cross-shaped valve have the following types of packaging:

a) in racks of 100 pcs per shrink film (in case the tubes are supplied with a capillary, a bag with 100 capillaries is attached externally to the tripod with adhesive tape). Tubes and capillaries are sterile. Transport packaging may contain 20, 25 or 30 tripods. Test tubes are delivered together with the application instruction (1 instruction for group packaging).

b) individually wrapped in PE film: mini-tube - 1 pc., capillary - 1 pc. (mini tube sterile packed with capillary). Individual packages are placed in a cardboard box containing 40 or 50 individual packages along with instructions for use (1 instruction for group packaging).

Transport packaging may contain 30 carton packs of 40 individual packs, or 24 carton packs of 50 individual packs.

Thus, a total of 1,200 individual packages fit in the transport packaging.

2. Disposable non-vacuum tubes for capillary blood with a screw cap have the following packaging:

- in tripods of 100 pcs per shrink film. Tubes are sterile.

Transport packaging may contain 20, 25 or 30 tripods.

Test tubes are delivered together with the application instruction (1 instruction for group packaging).

OPERATION, STORAGE, TRANSPORTATION CONDITIONS. SHELF LIFE

Keep non-vacuum tubes for capillary blood without blood at a temperature of 4-25 °C in a dry ventilated room. Keep away from harmful or poisonous substances.

Attention! Avoid direct sunlight. Exceeding the recommended storage temperature may lead to a loss of quality of the mini tube (for example, evaporation or staining of fillers (additives)).

Storage

Temperature	4-25 °C
Atmospheric pressure	86-106kPa
Relative humidity	Not more 75%

Operation

Temperature	15-30° C
Atmospheric pressure	86-106kPa
Relative humidity	35-75%

It is allowed to store samples in test tubes until testing for no more than 6 hours at a temperature of +18 to + 25°C or no more than 24 hours if stored from +2 to + 8°C and relative humidity from 30 to 75%, atmospheric pressure from 86 to 106 kPa

Transportation of products without blood should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the temperature from -25 to +25 ° C. During transportation to provide protection from the sun, rain and snow, do not mix with harmful or toxic substances.

Temperature	From -25 to +25 °C.
Air pressure	86-106kPa
Relative humidity	Not more 75%

For transportation of blood samples, it is necessary to use thermal containers that provide a temperature regime of +2 to + 8°C. Blood samples with known infectious material are placed in an additional secondary container. During transportation, the tubes with blood should be tightly closed and installed so as to prevent them from tilting.

Shelf life - 24 months

WARRANTY OBLIGATION

All questions related to the operation and quality of the product can be addressed to the Authorized representative of the manufacturer in the territory of the Russian Federation. The manufacturer guarantees that the product meets the requirements of the standards and documentation of the manufacturer under the conditions of transportation, storage, operation. The manufacturer shall not be liable for any damage or damage caused by the use of the product until proven guilty.

The manufacturer shall not be liable for any damage or damage caused by the use of the product after the expiry date indicated on the packaging.

For questions about the quality of a medical product, manufactured by HEBEI XINLE SCI&TECH CO.,LTD, China, please contact:

REAL CARE LLC

Perovskaya st. 1 building 10 floor 2 app VII room 3 Moscow Russia 111524

Telephone: +7(495)234-57-31

E-mail:lab@int-med.ru

POTENTIAL CONSUMER

User - Employee of a medical laboratory (doctor CDL or medical laboratory technician), medical staff (nurse) trained in taking biomaterial..

Patients – all ages

THE MAINTENANCE AND REPAIR OF MEDICAL DEVICES

Not applicable.

DISPOSAL OR DESTRUCTION OF MEDICAL DEVICES.

Used test tubes are potentially hazardous medical waste to the spread of infectious diseases transmitted through blood and are class B medical waste — epidemiologically hazardous waste. Disinfection and disposal of used tubes should be carried out in accordance with the requirements of sanitary regulations of SanPin 2.1.3684-21.

Unused tubes, including expired ones, should be disposed of as Class A waste and can be disposed of as household waste.

Adhere to national regulations when disposing or recycling the product, its components and its packaging.

A blood sample is a hidden source of contamination. Therefore, all components of the equipment, consumables and waste substances that come into contact with the sample should be used in accordance with biological safety rules.

Manufacturer:

Hebei Xinle Sci&Tech Co., Ltd China

No. 2, Xingye Street, Xinle City, 050700 Shijiazhuang City, Hebei Province People's Republic of China

+86 311 8858 1771

hbxinle001@hbxinle.net

Manufacturing site:

Hebei Xinle Sci&Tech Co., Ltd, Китай

No. 2, Xingye Street, Xinle City, 050700 Shijiazhuang City, Hebei Province People's Republic of China

AUTHORIZED REPRESENTATIVE IN RUSSIAN FEDERATION:

REAL CARE LLC

Perovskaya st. 1 building 10 floor 2 app VII room 3 Moscow Russia 111524

Telephone: +7(495)234-57-31

E-mail: lab@int-med.ru

APPENDIX

1. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
2. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
3. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
4. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
5. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
6. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
7. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
8. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
9. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
10. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
11. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
12. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
13. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with sodium heparin (NH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
14. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

15. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
16. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
17. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
18. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
19. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
20. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
21. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
22. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
23. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
24. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
25. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
26. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
27. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
28. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with lithium heparin (LH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

29. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
30. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
31. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
32. Mini tube non-vacuum for capillary blood (sterile, with lithium heparin (LH) (Green), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
33. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
34. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
35. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
36. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) . without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
37. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
38. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
39. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
40. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
41. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with rubber cap and cross-shaped valve,

- volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
42. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 43. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
 44. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 45. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 46. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
 47. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
 48. Mini tube non-vacuum for capillary blood (sterile, with sodium heparin (NH) and separating Gel (Green) , without capillary, with screw-down cover, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

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
证明书地址 Website for verifying the certificate: <http://www.rzceplt.com/validate.html>

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Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac))
 Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac))
 Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac))
 Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac))
 (see Appendix)

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In MicPlus tubes, a PP mini tube is placed inside the PET tube. The tube is made of a material that allows you to clearly see the contents of the container when visually inspected. A variety of the products allows you to pick up the suitable tube for specific laboratory research or analysis.

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The inner walls of the tubes and standalone capillaries are coated with micropieces of the anticoagulant EDTA K2 (2-substituted potassium salt of ethylenediaminetetraacetic acid) or EDTA K3 (3-substituted potassium salt of ethylenediaminetetraacetic acid). EDTA prevents blood from clotting by binding calcium ions and does not affect hematological parameters. A blood smear should be done within 3 hours after taking a sample, the rest of the tests are carried out for 6 hours at room temperature or for 24 hours at 0-2 ° C. Erythrocytes, leukocytes and platelets in a blood sample with an EDTA anticoagulant are stable for up to 24 hours.

Attention! If transportation of samples is required, tubes with EDTA K2 should be used, since EDTA K3 shows slightly lower properties as an anticoagulant and much more affects cell morphology during prolonged storage (more than 4 hours).

in a tripod at room temperature for at least 1 hour to form a strong gel barrier.

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Purpose: Test tubes are used for hematological, immunochemical, virological and other qualitative and quantitative studies of whole blood.

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1. All tubes for capillary blood collection of the required volumes and with the necessary fillers (additives).

- Capillaries (using tubes with a cross-shaped valve).
- 1. Lancets or scarifiers for puncture.
- 4. Personal protective equipment against contact with blood (disposable gloves, overalls, mask).
- 5. Wipes for cleaning the puncture site.
- 6. Dry sterile gauze or dry cotton swabs.
- 7. Container for recycling or its equivalent.

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1. Disinfect hands and wear disposable gloves.

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OR

Remove the capillary from the tube and insert into the tube of the second type through the cross-shaped valve on the cap (plug). If there is no capillary, open the cap to collect blood by gravity.

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Carefully remove the capillary from the minitube with a rubber cap and a cruciform valve (the valve of the stopper closes automatically). The used capillary should be immediately disposed of in a special container or its equivalent.

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Centrifuging

Make sure that the mini test tube caps do not rest on the walls of the centrifuge cap, otherwise the cover may come off. See the table below for sample preparation conditions for analysis. The recommended temperature is 15-25 ° C.

eparation of samples for analysis

Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
EDTA K2	24 hours	8-10 times	1300 g	10
EDTA K3	24 hours	8-10 times	1300 g	10

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Name	Volume	Dimensions	Additive	Cap (plug) color	Capillary	Cap type
MINI tube with EDTA K2 (lilac)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	K2EDTA	lilac	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
MINI tube with EDTA K3 (lilac)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	K3EDTA	lilac	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
Tube MicPlus with EDTA K2 (lilac)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	13x75 mm	K2EDTA	lilac	With a capillary and without a capillary	Rubber cap with cross-shaped valve
Tube MicPlus with EDTA K3 (lilac)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	13x75 mm	K3EDTA	lilac	With a capillary and without a capillary	Rubber cap with cross-shaped valve

MAIN RAW MATERIALS

Tube and capillary PP or PP+PET (MicPlus)
Rubber cap and cross-shaped valve Butyl rubber
Screw-down cover PE

STERILIZATION

The products are sterile (tubes and capillaries).
The individually wrapped test tube is sterilely packed together with the capillary.
The tube in a rack and shrink-wrapped is sterile.
The capillaries in the bag, which is attached to the outside of the rack, are also sterile.
Sterilization method - radiation sterilization.

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b) individually wrapped in PE film: mini-tube - 1 pc., capillary - 1 pc. (mini tube sterile packed with capillary). Individual packages are placed in a cardboard box containing 40 or 50 individual packages along with instructions for use (1 instruction for group packaging).
Transport packaging may contain 30 carton packs of 40 individual packs, or 24 carton packs of 50 individual packs.

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E-mail:lab@int-med.ru

POTENTIAL CONSUMER

User - Employee of a medical laboratory (doctor CDL or medical laboratory technician), medical staff (nurse) trained in taking biomaterial.

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E-mail:lab@int-med.ru

APPENDIX

1. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
2. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
3. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
4. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
5. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
6. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
7. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
8. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
9. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
10. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
11. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
12. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
13. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
14. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)

15. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
16. Mini tube non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
17. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
18. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
19. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
20. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
21. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
22. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
23. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
24. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
25. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
26. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
27. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
28. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

29. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
30. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
31. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
32. Mini tube non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
33. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
34. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
35. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
36. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
37. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
38. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
39. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
40. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K2 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
41. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, объем volume of capillary blood sample 0,2 ml. Sizes 13 mm x 75 mm.
42. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)

43. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
44. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs in rack)
45. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
46. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
47. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)
48. Tube MicPlus non-vacuum for capillary blood (sterile, with EDTA K3 (lilac), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 13 mm x 75 mm, quantity in pack 100 pcs, in rack)

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THIS IS TO CERTIFY THAT: the seal of HEBEI XINLE SCI & TECH CO.,LTD. on the annexed INSTRUCTION FOR USE is genuine.

China Council for the Promotion
of International Trade

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
证明书查询地址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red),
 Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), (see Appendix)
 RC № RZN XXX from XX.XX.XXXX

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Name of the tube ¹	Cap (plug) color
MINI tube with Clot-activator (Silica)	Red
MINI tube with Clot-activator (Silica) and separating Gel	yellow

The inner walls of the tubes are coated with microparticles of silicon oxide, which activate the process of the formation of a clot by gently inverting the tube. The recommended clotting time is not less than 10- 30 minutes.

Contain a separating gel at the bottom of the tube. The density of this material is lower than that of a blood cell, but higher than that of serum. When centrifuged, the gel moves to the border between the blood cell and the serum, where it forms a stable barrier. Due to this, the serum

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can be sampled directly from the mini tube. Serum remains stable for up to 48 hours under recommended storage conditions in gel tubes.

Attention! Gel tubes cannot be transported without prior centrifugation! Sharp movements can destroy the gel and lead to the ingress of its particles into the sample, so that the results of the analysis will be erroneous. If transportation is required, then after centrifugation, it is necessary to stand the test tubes with gel in a tripod at room temperature for at least 1 hour to form a strong gel barrier.

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Purpose: for obtaining and research of serum. Test tubes are used for biochemical, serological, immunological and microbiological studies, as well as for drug monitoring

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1. All tubes for capillary blood collection of the required volumes and with the necessary fillers (additives).
2. Capillaries (using tubes with a cross-shaped valve).
3. Lancets or scarifiers for puncture.
4. Personal protective equipment against contact with blood (disposable gloves, overalls, mask).
5. Alcohol wipes for cleaning the puncture site.
6. Dry sterile gauze or dry cotton swabs.

- Container for recycling or its equivalent.

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1. Disinfect hands and wear disposable gloves.
2. Prepare a mini test tube:
Open the cap of the screw cap tube designed for gravity collection of blood by turning it counterclockwise.
OR
Remove the capillary from the tube and insert into the tube of the second type through the cross-shaped valve on the cap (plug). If there is no capillary, open the cap to collect blood by gravity.
3. Disinfect the puncture site. Prepare the automatic lancet or scarifier according to the manufacturer's instructions.
Attention! Do not use cotton swabs for disinfection, as their fibers can get into the sample and lead to incorrect analysis results (for example, clog the cell count sensor and hemoglobin chamber).
4. Pinch the patient's finger between the thumb, middle and index fingers. OR
For newborns: gently, but firmly take the child by the heel.
5. Make a puncture, remove the first protruding drop.
6. Gently press down on the surrounding tissue so that a drop of blood appears. Without touching the puncture site, collect the dripping blood:
By gravity into a test tube with a cap that turns with the help of a tray ledge. OR
Through the capillary into the mini-test tube with a rubber cap and a cruciform valve (the minitube should be held horizontally). Wait until the capillary is filled with blood and turn the tube vertically so that the blood can freely enter it from the capillary.
Attention! Watch out for the filling level of the minitube. Compliance with the required level of filling allows you to achieve the correct ratio of blood and filler (supplements).
7. After filling the tube:
Close the screw cap of the mini test tube tightly. OR
Carefully remove the capillary from the minitube with a rubber cap and a cross-shaped valve (the valve of the stopper closes automatically). The used capillary should be immediately disposed of in a special container or its equivalent.
Attention! If you collect blood in several mini tubes, a new capillary should be used for each.
8. Gently invert the mini-tube 5-7 times to ensure thorough mixing of the filler (additive) and blood. Tap lightly on the bottom of the microfuge to ensure that the filler (additive) and blood are well mixed.
Attention! Do not shake the tubes, as this may cause churning or hemolysis. Insufficient or slow inverting the serum tubes can slow down clot formation. In tubes with anticoagulants, improper mixing can lead to gluing of platelets, clot formation and incorrect results of the analysis.
9. When all necessary samples have been collected, press a dry sterile swab into the puncture site and wait until the bleeding stops. (If blood was taken from the heel, lift the baby's leg up.)
10. After the blood is drawn, the traces of blood may remain in the recess of the mini-vial cap with the rubber cap and the cross-shaped valve. Take precautions to avoid contact with it.

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Additive	Maximum time after blood sampling	How many times to flip	Recommended relative centrifugal acceleration	Recommended centrifugation time, min
Clot activator	3 hours	5-6 times	1500 g	10
Clot activator+gel	3 hours before centrifugation and 48 hours after centrifugation	5-6 times	1800 g	10

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GUARANTEED BY MANUFACTURER VALUES OF MAIN PARAMETERS, CHARACTERISTICS (FEATURES) OF MEDICAL DEVICE.

Values of the parameters, characteristics (features) of the medical device are presented in the table below. Where applicable, the permissible deviations from the nominal characteristics are $\pm 10\%$, unless otherwise agreed.

Name	Volume	Dimensions	Additive	Cap (plug) color	Capillary	Cap
MINI tube with Clot-activator (Silica) (RED)	0,2 ml, 0,25 ml, 0,5 ml, 1 ml	10x45 mm	Clot-activator (Silica)	Red	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover
MINI tube with Clot-activator (Silica) and separating Gel (Yellow)	0,2 ml, 0,25 ml, 0,5 ml, 0.75 ml	10x45 mm	Clot-activator (Silica) and Gel	Yellow	Rubber cap with cross-shaped valve - With a capillary and without a capillary Screw-down cover – without capillary	Rubber cap with cross-shaped valve/ screw-down cover

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E-mail:lab@int-med.ru

APPENDIX

1. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
2. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
3. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
4. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs. in rack
5. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
6. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
7. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
8. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
9. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in ack)
10. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
11. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
12. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
13. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
14. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)

15. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
16. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) (red), without capillary, with screw-down cover, volume of capillary blood sample 1 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
17. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm. , quantity in pack 100 pcs in rack)
18. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
19. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
20. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with screw-down cover, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
21. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
22. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
23. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
24. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with screw-down cover, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
25. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
26. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

27. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
28. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with screw-down cover, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
29. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs in rack)
30. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
31. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
32. Mini tube non-vacuum for capillary blood (sterile, with Clot-activator (Silica) and separating Gel (Yellow), without capillary, with screw-down cover, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)

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

证明网址 Website for verifying the certificate: <http://www.rzccpit.com/validate.html>

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Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) (see Appendix)

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All tubes and capillary are sterile, equipped with a convenient label with a field for recording, which helps to easily identify the tube. Each tube has a calibration mark to accurately determine the level of filling (filling line) located on the label or on the tube itself, depending on the intended volume. Compliance with the required level of filling allows to achieve the correct ratio of blood and filler (additive).

Name of the tube ¹	Cap (plug) color
MINI tube (Si-plus) with sodium citrate 4: 1	Blue

Si plus tubes contain a 3.8% buffered solution of trisodium citrate (0.129 mol / L). The blood reagent ratio is 4: 1. Test tubes are used for syphilis screening.

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Purpose: for serological tests for syphilis of capillary blood.

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1. Disinfect hands and wear disposable gloves.
2. Prepare a mini test tube:
Remove the capillary from the tube and insert into the tube of the second type through the cross-shaped valve on the cap (plug). If there is no capillary, open the cap to collect blood by gravity.
3. Disinfect the puncture site. Prepare the automatic lancet or scarifier according to the manufacturer's instructions.
Attention! Do not use cotton swabs for disinfection, as their fibers can get into the sample and lead to incorrect analysis results (for example, clog the cell count sensor and hemoglobin chamber).
4. Pinch the patient's finger between the thumb, middle and index fingers. OR
For newborns: gently, but firmly take the child by the heel.
5. Make a puncture, remove the first protruding drop.

5. Gently press down on the surrounding tissue so that a drop of blood appears. Without touching the puncture site, collect the dripping blood:
By gravity into a test tube with a cap that turns with the help of a tray ledge. OR
Through the capillary into the mini-test tube with a rubber cap and a cross-shaped valve (the minitube should be held horizontally). Wait until the capillary is filled with blood and turn the tube vertically so that the blood can freely enter it from the capillary.
Attention! Watch out for the filling level of the minitube. Compliance with the required level of filling allows you to achieve the correct ratio of blood and filler (supplements).
7. After filling the tube:
Carefully remove the capillary from the minitube with a rubber cap and a cruciform valve (the valve of the stopper closes automatically). The used capillary should be immediately disposed of in a special container or its equivalent.
Attention! If you collect blood in several mini tubes, a new capillary should be used for each.
8. Gently invert the mini-tube 5-7 times to ensure thorough mixing of the filler (additive) and blood. Tap lightly on the bottom of the microfuge to ensure that the filler (additive) and blood are well mixed.
Attention! Do not shake the tubes, as this may cause churning or hemolysis. Insufficient or slow inverting the serum tubes can slow down clot formation. In tubes with anticoagulants, improper mixing can lead to gluing of platelets, clot formation and incorrect results of the analysis.
9. When all necessary samples have been collected, press a dry sterile swab into the puncture site and wait until the bleeding stops. (If blood was taken from the heel, lift the baby's leg up.)
10. After the blood is drawn, the traces of blood may remain in the recess of the mini-vial cap with the rubber cap and the cross-shaped valve. Take precautions to avoid contact with it.

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Transportation of products without blood should be carried out by all types of covered transport in compliance with the conditions and requirements established for this type of transport, at the temperature from -25 to +25 °C. During transportation to provide protection from the sun, rain and snow, do not mix with harmful or toxic substances.

Temperature	From -25 to +25 °C.
Air pressure	86-106kPa

Relative humidity	Not more 75%
-------------------	--------------

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E-mail:lab@int-med.ru

APPENDIX

1. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
2. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
3. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,2 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
4. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
5. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
6. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary with rubber cap and cross-shaped valve, volume of capillary blood sample 0,25 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
7. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,375 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
8. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,375 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
9. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,375 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
10. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
11. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
12. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,5 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)
13. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , without capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
14. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, quantity in pack 100 pcs, in rack)
15. Mini tube Si-plus non-vacuum for capillary blood (sterile, with sodium citrate 4: 1 (Blue) , with capillary, with rubber cap and cross-shaped valve, volume of capillary blood sample 0,75 ml, sizes 10 mm x 45 mm, individually wrapped, 40 or 50 individual packs in carton)

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