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Sales Management Russia / CIS Region Europe
_____ /Martin Koch/

Regulatory Affairs Manager Eastern Europe
_____ /Olga Timofejev/

Dräger

Dräger Safety AG & Co. KGaA
Revalstrasse 1
23560 Lübeck, Germany
www.draeger.com

_____ 2016

User Manual

Breath Alcohol Analyzer Alcotest 5510
with Accessories

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1. INTRODUCTION

Prior to working with the Alcotest 5510 analyzer of ethanol va in the expiratory air (hereafter referred to as the "analyzer"), please read and follow the instru of this operating manual.

This operating manual contains complete information about the alyzer's functioning principles, operating procedures, capabilities, specifications, frequency nt of maintenance.

2. DESCRIPTION AND OPERATION

2.1 Purpose

The Alcotest 5510 analyzer with accessories is used for det ining an ethanol vapor's mass concentration in the sampled expiratory air in medical facilities and for personal use.

2.2 Specifications

Analyzer Alcotest 5510 is produced and discharged witl tol nces for dimensions and weight, according to the specifications given in Table 1.

Table 1

Parameter	Specifications Alcotest 5510
Basic dimensions and parameters:	
Instrument type	Portable digital breath alcohol analyzer
Sensor type	Electrochemical
Display	Graphic, monochrome
Display results	Numeric, three digits
Additional display	Audible alarms
Product weight, g	195
Overall dimensions, H × W × D, mm	142 × 75 × 36
Duration of preparations for the first test	No more than 6 seconds, with cut-off automatic control of the absence of ethanol in the ambient air; No more than 20 seconds with enabled automatic control of the absence of ethanol in the ambient air;
The duration of the preparation for re-test	No more than 40 seconds after the analysis of the gas mixture with the mass concentration of ethanol was 0.50 mg / l
The minimum amount of exhaled air samples	At least 1.2 liters.
The test mode	1. Measuring mode— with expiration through a disposable mouthpiece and automatic sampling 2. Screening mode - expiration without mouthpieces with manual sampling 3. Manual mode— with expiration through a disposable mouthpiece and manual sampling
Presentation of results	1. Measuring mode—a three-digit number 2. Screening mode-- a text message, "Alcohol" or "No

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

	alcohol" 3. Manual mode—a thre number
Blowing time	4 s
Memory, tests	10 tests, without specify date and time of the test
Ability to connect to the printer to print the results on paper	No
Ability to connect to the computer	No
Mode of operation unit	Intermittent
The sound level of audio signals	At least 40 dB
The term electrochemical sensor life	At least 2 years
Mean lifetime	6 years
Mean time between failures	16000 hours
Identification data (codes)	
Software identification name	8323547
Software revision number (identification number)	Rev.1.10
Digital identifier of software	0x6202 (CRC16)
Electrical specifications:	
Instrument power supply	Alkaline AA batteries (or ble NiMH batteries
Number of batteries or accumulators	2 pcs.
Operating voltage	3.0 V DC (batteries) or 2.4 V DC (accumulators)
The number of tests with fully charged battery, not less	1500
Accumulator charging	Using a self-contained ch device
Requirements for environmental conditions	
Conditions of use:	
Place of use	Indoor and outdoor use
Ambient temperature	from – 5°C to + 50°C
Relative humidity	from 10% to 100% (with ndensation)
Atmospheric pressure	from 84 kPa to 106.7 kPa
Conditions of transportation and storage	
Ambient temperature	from – 40°C to + 70°C
Relative humidity	from 10% to 100% (without ndensation)
Atmospheric pressure	from 60 kPa to 140 kPa
Special requirements for conditions of transportation and storage	Keep the instrument away fr n alcohol-containing liquids
Reliability requirements	
Climatic performance	UHL 4.2
The degree of water protection	IP54
Resistance depending on sensed mechanical impacts	For Group 3

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Safety requirements	
The extent of the potential risk of application	Class 1
Degree of protection against electric shock	The product with the [illegible] pply, Type B
The degree of security against flammable anesthetics or mixtures of oxygen	It is not suitable for [illegible] anesthetic mixtures containing air, oxygen [illegible] xide
The metrological characteristics	
Measurement range	0.00–2.00 mg/l
Indication range	0.00–2.50 mg/l
Resolution of minimal scale discharge	0.01 mg/l
The limits of maximum permissible errors	According to Table 2
The limits of maximum permissible errors, due to changes in external factors affecting the workers within the operating conditions relative to normal conditions	According to Table 3
Additional error of the presence not measured components	Missing
Limit of acceptable standard deviation	No more than 1/3 as a frac[illegible] of the maximum permissible error analyzers at a temper[illegible] 5) ^{0}C
Maximum permissible variation of readings	No more than 0.5 as a fracti[illegible] maximum permissible error analyzers at a temperat[illegible] (20) ± 5) ^{0}C
The interval of the analyzer readings without correction	Not less than 12 months

The limits of maximum permissible errors of th[illegible]lyzer of ethanol vapor in the breath Alcotest 5510 depending on the r[illegible]ge of concentrations measured at ambient temperature [illegible]) ± 5) ^{0}C.

Table 2

The measuring range of the mass concentration of ethanol in mg / l	The maximum permissible [illegible]s at temperature (20 ± 5) ^{0}C	
	absolute	relative
0-0,20	±0,02 mg / l	-
over 0,20-2,00	-	±10%
Remarks: 1) Measurement on "mg / l" mass concentration of ethanol in th[illegible]isplay unit displays a «mg / l». 2) When testing the analyzer using a gas mixture the composit[illegible] of ethanol [illegible] nitrogen cylinders under pressure using the coefficient allocation indications K[illegible]. [illegible]l to 1.05.		

The limits of maximum permissible errors of t[illegible] analyzer of ethanol vapor in the breath Alcotest 5510, depending on the ambient temperature

Table 3


Environment	The measuring range	The maximum permissible errors [1]

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

temperature	of the mass concentration of ethanol in mg / l		relative
From minus 5.0 °C to 15.0 °C incl.	0-0,33		-
	over 0,33-2,00		±15%
St. 15.0 °C to 25.0 °C incl.	0-0,20	±0.	-
	over 0,20-2,00		±10% 2)
St. 25,0 °C to 50,0 °C incl.	0-0,33		-
	over 0,33-2,00		±15%
1) The table shows the maximum permissible error anal			nditions provided in
the "Description of the type of measuring instruments", p			
2) As shown in Table 2.			

2.3. Product Components

Figure 1 shows the opened case with the analyzer and acc as a standard complete set.

st analyzer
posable mouthpieces
V size piles or batteries
nd strap
for transportation and storage
erating manual

Figure 1—Case with the analyzer and accessories

1. Alcotest analyzer

The Alcotest 5510 analyzer housing is made of ABS pla me: PA-77NB) resistant to high humidity, temperature, as well as static electricity. It contai electronic board with a display, electronic circuit elements, and contact pads for batteries; trochemical ethanol sensor; a pressure sensor; an electromagnetic pump for air sampling.
Dimensions (H×W×D): 142×75×36 mm. Weight (with electr lls): 195 g.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

2. Plastic case

The case is made of polypropylene (grade: 1500J). A compartment is provided for each **element** inside it to prevent them from impacts and other mechanical damage. **Dimensions** (H×W×D): 233×178×55 mm. Weight: 292 g.

3. Hand strap


The hand strap is designed to prevent the instrument falling so as to avoid any negative consequences.
Material: 100% polyester (grade: 600D Rip Stop). Length: . Weight: 1 g.

4. Operating manual (hard copy, in Russian).

Accessories:


1. Mouthpieces

The disposable mouthpieces are designed for an exhale with an accurate result. For this purpose, they are attached directly to the analyzer housing.
Material: food-grade polystyrene (grade: 158K). Dimensions (H×W×D): 110×24×15 mm. Weight: 4 g.

2. Mouthpieces non-return valve

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The mouthpieces with non-return valve are designed for exhalation test with an accurate result. For this purpose, they are attached directly to the analyzer housing and prevent the person being tested from inhaling air.

Material: food-grade polystyrene (grade: 158K). Dimensions (H×W×D): 110×24×15 mm. Weight: 4 g.


3. Alkaline batteries or rechargeable batteries

The alkaline batteries (rated voltage: 1.5 V) or rechargeable batteries (rated voltage: 1.2 V) are intended for self-contained power supply to the analyzer.

Dimensions (H × D): 50 × 14 mm. Weight: 23 g.

4. Leather bag

The leather bag is designed to attach the analyzer to a trouser belt, as well as to protect the device against the environment.

Material: 100% natural leather (splitleather). Dimensions (H×W×D): 175×70×25 mm. Weight: 106 g.

5. Sticker

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Sticker are designed to be attached to the rear panel of the analyzer.
Dimensions (H×W): 80×60 mm. Weight: 2 g.

6. Charger

The charger is designed to recharge the AA accumulators from 220 V mains.
Housing material: ABS plastic (grade: PA-77NB). Dimensions (H × W × D): 110×90×46 mm. Weight: 199 g.

The analyzer's outward appearance is shown in Figure 2.

Figure 2—The analyzer's outward appearance

There are three control buttons and a graphical monochrome display on the front panel of analyzer.

The back side of analyzer is provided with a battery compartment cover.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

2.4. Design and Operation

The Alcotest 5510 analyzer is a portable electronic device that quickly and accurately determines an ethanol vapor's mass concentration in the expiratory air. For operator's convenience, messages about all preparation and testing stages are output on a full-text display in Russian and followed by an audible signal.


The analyzer Alcotest 5510 is equipped with an original electrochemical ethanol sensor which provides reliable readings and high speed of operation. Ethanol vapors in exhaled air get in contact with the sensor surface and cause a reaction, which results in generating an electrical signal by the sensor. The chemical formula of the process, describing the operating principle of an electrochemical sensor for alcohol is as follows:

$C_2H_5OH + O_2 \rightarrow C_2H_4O_2 + H_2O$

Ethanol electrochemical sensor consists of two active surfaces platinum (electrode), which is deposited on a specially selected catalyst in the presence of which weaken bonds in ethanol. Between the electrodes is a porous substrate, impregnated with electrolyte. When the air blowing containing ethanol vapor, through the body sensory cells from ethanol molecules near the surface of the active catalyst come off the positively charged hydrogen ions H + and free electrons are formed. H + ions migrate through the electrolyte falling on the second surface, and wherein the water formed by contacting with oxygen. Thus, on one surface formed due to excess negative charge of the electron, on the other - positive. Being included in the electric circuit, the sensor circuit generates a current whose strength at the peak is proportional to the amount formed on the surfaces of the free charge, and hence the number of ethanol molecules reacted near the active catalyst surface.

The shape change of the sensor signal versus time is shown in Fig. 2

Figure 2. The shape changes of the electrochemical sensor signal

Sensor recovery occurs naturally as the decay of a chemical reaction, and it is the more, the higher the concentration of ethanol in the air sample, was admitted to the sensor. Built-in microprocessor that controls the entire operation of the analyzer converts the output of the sensor signals in the analyzer are displayed on the graphic display.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The external control of the analyzer is provided using buttons on the front panel. In order to facilitate operation, the instrument also has an audible signal system. Power supply is provided by 2 AA batteries or accumulators, and their charge is sufficient for up to 1,500 tests.

The analyzer Alcotest 5510 has three modes of operation in accordance with different types of blowing:

1. Measuring mode: direct blowing through a disposable mouthpiece within a specified time. The instrument automatically takes an air sample for analysis upon reaching the specified criteria of exhalation—sample volume and exhalation duration. The result of such measurement is presented on a display as a digital value with resolution to two decimal places.

2. Screening mode: an exhaled air sample is taken without mouthpieces. In this case, the result is an estimate and is presented in the form of "Alcohol" or "No alcohol" messages.


3. Manual mode: air is exhaled through a mouthpiece, but a sample is taken manually by an operator without waiting until the end of exhalation. The result of such measurement is presented on a display in digitized form with the comment "Manual."

When using the instrument, you should have in mind that the accuracy of measurements is only ensured in the measuring mode.


The analyzer Alcotest 5510 has a memory for storing results and a menu system allowing the operator to view the stored results of last 10 tests in order and date of last calibration.

Each time the analyzer performs a self-test and monitoring of environmental parameters, when it detects an internal failure or deviation of any parameter is outside the permitted range, the instrument gives the display a warning about the detected fault or error. Until the cause of a malfunction or error can not work with the analyzer.

2.4.1. Outward Appearance and Description of the Main Elements of Analyzer

1—Display
2—"Up" menu button
3—"OK" button
4—"Down / Menu" menu button
5—Connector
6—Red / green / yellow indicator
7—Mouthpiece holder

The back side of analyzer is provided with a space for "stickers with a quick start guide".

After acquisition of the analyzer:

- Detach the sticker from the foil base.
- Place the sticker on its intended area on the back side of housing.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The symbols below the latch of battery compartment cover mean:
ATTENTION! Read the operating manual!

2.4.2. Display

In the measuring mode, the display shows a measurement result, as well as information about running processes, and relevant menu items are highlighted at the menu access area of the analyzer. The display is back-lighted.

2.4.3. Buttons

The ▼"DOWN" and ▲"UP" buttons are used for entering values and selecting menu functions.

The "OK" button confirms entered parameters or selected functions and is used for switching the analyzer and forced air sampling (in a manual sampling mode) on/off.

2.4.4. Indication

The three-color indicator complements readings on the display.

2.4.5. Menu Language

The default language is Russian.

2.4.6. Mouthpiece Installation Method

A nipple with an inlet opening is located in the recess of mouthpiece holder. For added convenience, the design of analyzer allows for both the left and the right orientation of mouthpieces with respect to the front panel.

2.4.7. Disposable Mouthpieces

The individually packaged changeable disposable plastic mouthpieces are used for breathing out into the analyzer. During breathing out through the mouthpiece, the analyzer automatically controls the duration and air flow rate in accordance with the factory settings. A sample of air is taken at the end of breathing-out that makes it possible to measure an ethanol content in the air from the depth of lungs.

2.4.8. Internal Memory

The analyzer has an internal memory that stores the results of the last 10 measurements with a test sequence number. The memory is accessed through the user menu.

2.4.9. Operation Modes

The analyzer has two air sampling modes—the automatic and manual ones. In order to take samples of air in the automatic and manual modes, the changeable disposable plastic mouthpieces are used, and measurement results are output on the display in digitalized form.

Moreover, the analyzer allows for air sampling in a passive mode, without a mouthpiece, in which case an assessment of ethanol presence in the air is made, and the display shows messages "Alcohol" or "No alcohol."

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

2.4.10. Software

The analyzer has a firmware. The analyzer firmware is developed by the manufacturer specifically to solve the problem of measuring an ethanol vapor's mass concentration in the expiratory air, as well as to show measurement results on the display and to store measured data. A firmware identification of the analyzer is made in the service menu by outputting a version on the analyzer display.

3. INTENDED USE

The analyzer Alcotest 5510 with accessories is applied for measuring the weight concentration of ethanol vapors in a collected sample of exhaled air at health care facilities and for personal use.

The analyzer Alcotest 5510 is intended for professional use by the physicians and paramedics who perform medical examination for alcohol intoxication, as well as by health care workers in rooms for pre-trip medical examination of drivers. It can be used as a work tool by emergency medical service providers.

The analyzer Alcotest 5510 is a hand-held portable instrument with self-contained power supply from two AA batteries or accumulators. It is equipped with a highly-selective and durable electrochemical ethanol sensor and provides a quick start and short analysis time after testing. The obtained results can be stored in a memory and easily read on a graphic display.

The analyzer Alcotest 5510 has a control menu with a two-level structure: a first-level menu for an operator and a second-level menu protected by a PIN code for service engineers.

3.1. Indications

Indications for use the Alcotest 5510 analyzer is the need to check for alcohol in the human body, particularly in case of the following external manifestations:

1. Alcohol vapors from mouth
2. Posture instability
3. Speech disturbance
4. Apparent fingers tremor
5. Sudden change in color of the facial skin
6. Behavior inappropriate to situation

3.2. Contraindications

Contraindications to use the Alcotest 5510 analyzer do not exist.

3.3. Possible Side Effects

Side effects during the use of the Alcotest 5510 analyzer do not exist.

3.4. Preparation of the Product for Use

Operating Limitations

- Ambient air must be free of alcohol and solvents (disinfectants) vapors.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

- Do not expose the analyzer sensor to the influence of gas mixtures with a high ethanol content, for example, by rinsing the mouth cavity with alcohol immediately before a measurement. This reduces the service life of sensor.
- It is not allowed to store and use the analyzer in premises wherein alcohol-containing substances are stored or work is carried out with them.
- Do not bring the analyzer close to antennas of mobile phones and transmitting stations.
- Repairs of the analyzer should be performed by qualified specialists in service centers.
- A new disposable mouthpiece should be used for each examinee.

3.4.1. Visual Inspection of the Analyzer

- Check for a seal (a special sticker on the fixing screw of the analyzer's back panel, with the battery compartment cover removed) and operability of the controls.
- Make sure there is no mechanical damage.

3.4.2. Operating Conditions

Before using, keep the analyzer in the operating conditions for 1 hour if the storage conditions were not consistent with the operating environment.


During the operation must be maintained following parameters:


Temperature: - 5°C to + 50°C

Relative Humidity: 10 - 100%

Atmospheric pressure: 84 – 106,7 кPa

3.4.3. Batteries Installation

When preparing the analyzer for operation or if the display shows a low battery symbol, install two AA alkaline batteries (2×1.5V) or two AA rechargeable batteries (2×1.2V) into the analyzer:

- Carefully press the latch of battery compartment cover and remove the cover.
- Install the batteries while observing the polarity indicated on the housing of analyzer.
- Close the battery compartment cover.

Battery charging:

In order to recharge the AA size batteries, remove them from the analyzer and put them in a special battery charger. Follow the operation guidelines for the battery charger.

3.4.4. Checking the Battery Charge

A battery charge status can be determined by the battery symbol in the lower right corner of the display. The analyzer has four battery charge indication levels:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU


Full charge

⅔ of the total capacity of batteries

⅓ of the total capacity of batteries

The batteries are discharged—in this case, replace the piles with new ones or charge the batteries.


3.4.5. Fastening the Strap onto the Hand


To eliminate the risk of accidental fall of the analyzer, it is possible to attach the strap to the hand and put it on the wrist during operation:

- Insert the small loop into the housing hole.
- Thread the long loop into the small loop and tighten it.

3.4.6. Mouthpiece Installation

- Open the packaging. For hygienic reasons, do not remove the package from the section being in contact with the mouth until the mouthpiece is properly installed into the analyzer.
- Slide the mouthpiece diagonally (the riffled side upwards) along the guiding groove until it stops, and press the mouthpiece into the holder. A clicking sound will confirm the correct engagement.
- Completely remove the packaging.

The mouthpiece may be oriented to both the right and the left side.

3.5. Operation Procedure

When performing a test, observe the following rules:

a) An analyzable air sample must contain neither tobacco smoke particles, nor alcohol residues or alcohol-containing medicated preparations from the mouth cavity, nor sputum and saliva.

Therefore, prior to carrying out a test:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

– At least 2 minutes must have elapsed after smoking.
– Not less than 20 minutes must have elapsed after use of alcohol, alcohol-containing medicinal preparations, oral sprays, as well as food containing low concentrations of alcohol (fermented milk products, kvass, etc.).

ATTENTION! Rinsing the mouth cavity with water or soft drinks does not reduce the necessary waiting period!

b) Before a test, the examinee must breathe normally, not perform a hyperventilation (a rapid tidal breathing).

c) During a test, the examinee must provide the required flow rate and volume of expiratory air. The expiratory air flow must be constant (non-stop).

3.5.1. Attaching the Mouthpiece

Insert a new mouthpiece into the analyzer holder.

3.5.2. Turning on the Analyzer

Press and hold the "OK" button for about 1 second until the display shows the start screen.	

3.5.3 Preparing the Analyzer

After 2 seconds the analyzer will enter the self-testing mode, which takes about 5 seconds and is accompanied by the text **Самотестир (Self-testing)** on the display. Then the analyzer will automatically take a sample of ambient air (if the automatic control function of absence of ethanol in the ambient air is on). At the same time, a click is heard, a message **нуль тест (zero test)** is displayed a few seconds, and if alcohol is not detected in the ambient air, a message **ГОТОВ (READY)** is output on the display, and a short beep sounds. Simultaneously, a number of current test is output, and a green indicator flashes.

If alcohol is detected in the ambient air, a message **Обнаружен алкоголь. Прибор не готов! (Alcohol is detected. The device is not ready!)** is output on the display, and a red indicator is on. Wait for an automatic shutdown of analyzer and turn on the analyzer again.

If the analyzer is turned off forcedly, the problem will persist when restarting. In order to reset, it is necessary to remove the batteries from the analyzer for a few seconds, and then to install them back.

3.5.4. Performing a Test

• After a message **ГОТОВ (READY)** is output on the display, ask the examinee to breathe out into the wide entrance of mouthpiece. The breathing-out should be vigorous, without forcing and interruption. A correct breathing-out is designated by a continuous beep and a lighted green indicator.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

• During sampling, the display shows a message: **ДУЙТЕ (BREATHE OUT).**


After sampling, the green indicator goes out, and the continuous beep is switched off.

• The display shows: **Идет анализ (Analyzing).**

• After 5–25 seconds (depending on temperature and a measured concentration value), a measurement result of ethanol mass concentration in the expiratory air is output on the display as a message **X.XX mg/L**.


The units of measurement of the ethanol mass concentration "mg/L" are displayed as **mg/L**.

In addition, simultaneously with a displayed measurement result, a beep sounds: 1 time—low level, 2 times—middle level, 3 times—high level.


3.5.5. Errors during the Testing

• **If** a volume of supplied sample is below an allowable value (see 1.2.6), the display shows a message: **Недостаточный объем (Insufficient volume),** a short beep sounds, and a red indicator flashes.

To repeat the test, press the "OK" button, and the analyzer will be ready for a new test in 4 seconds.


• If a flow rate of supplied sample is below an allowable value (see 1.2.6), in case of an interruption of breathing-out or a sucking in at the end of sampling, the display shows a message **Выдох прерван (Breathing-out is interrupted),** a short beep sounds, and a red indicator flashes.

To repeat the test, press the "OK" button, and the analyzer will be ready for a new test in 4–30 seconds.

3.5.6. Manual Sampling Mode

If the examinee fails to make a breathing-out that triggers an automatic sampling, *i.e.,* he/she cannot provide the required expiratory flow rate and volume (see 1.2.6), it is possible to use a manual sampling mode.

To this end, press the "OK" button when the examinee breathes out, in which case the analyzer will perform a manual sampling. Perform a manual sampling at the very end of breathing-out which the examinee is capable of.

A message **Ручной (Manual)** appears below the displayed measurement result.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU


ATTENTION! During the manual sampling mode, breathing-out parameters are not provided. During the manual sampling mode, a measuring inaccuracy may exceed a maximum permissible error of analyzer.

3.5.7. Removing the Mouthpiece

- Press the mouthpiece from below, as shown in the figure.
- Take the mouthpiece out of the holder.
- Dispose of the mouthpiece in accordance with local regulations.

3.5.8. Preparation for the Next Testing

In order to prepare for the next testing, press the "OK" button. During cleaning of the sensor, the display counts down in seconds and shows a rotating icon.

Note: Turning the analyzer off and on at this stage does not reduce the preparation time of analyzer for operation.

3.5.9. Turning off the Analyzer

- Press and hold the "OK" button. The analyzer will turn off in 2 seconds.
- An automatic shutdown of analyzer takes place in 4 minutes after a stoppage of work.
- The analyzer automatically shuts off in case of a low voltage of piles / batteries power supply. At the same time, a battery symbol flashes on the analyzer display, and in order to continue, it is necessary to replace the power supply piles or recharge the batteries.

3.5.10. Operation Procedure in the Passive Sampling Mode

ATTENTION! The passive sampling mode can be used for a preliminary assessment of ethanol presence / absence in the air. This mode does not provide for metrological characteristics of the analyzer.

- Keep the analyzer in the analyzable air, without installing the mouthpiece.
- Briefly press the "OK" button. A click will be heard, then an analysis will start.
- An assessment result will appear on the display as a message:
 - a) **Алкоголя нет (No alcohol)**
 - b) **Алкоголь (Alcohol)**

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

• Press the "OK" button for performing the next testing.

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3.5.11. Operation Procedure in the Menu of Analyzer

3.5.11.1. Menu Symbols

A transition to a higher menu level
A menu item that contains a submenu
An open menu item
A selected function
There are menu items below visible items
There are menu items above visible items
An activated item
An activated screen part

3.5.11.2. Menu Access

Turn on the analyzer. Menu functions can be called when a message **ГОТОВ (READY)** is displayed.

• In order to access the user menu, press and release the button ▼.

• In order to access the service menu (it is intended only for specialists of service centers and verification officers), press and release the button ▼, then enter a PIN code.

3.5.11.3. Operation in the Menu

• Select a function using the buttons ▼ and ▲.

• The selected function is highlighted.

• Activate the selected function by pressing the "OK" button.

3.5.10.4. Exit from the Menu

• After selecting the line **Возврат в изм. режим (Back to the measuring mode)** and pressing the "OK" button, it is possible to switch the analyzer over to the measuring mode.

• If no buttons are pressed during 2 minutes, the analyzer automatically returns to the measuring mode.

3.5.10.5. User Menu Items

Последн. тестов (Of the last tests)	The results of the last 10 measurements are output on the display, data selection is carried out by the buttons ▼ and ▲.
Дата посл. калибр. (Date of the last calibration)	The date of the last readings adjustment is output on the display.

Returning from a selected user menu item is accomplished by the "OK" button.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

4. DISINFECTION

After each use, should wipe clean the outside of the Analyzer with a damp cloth soaked in soapy water. While taking care of the analyzer must not be used alcohol-based, abrasive cleaners and tools with dyes, as well as to carry out processing with hard brushes.
Mouthpieces are disposable and cannot be reused and processed.

5. MAINTENANCE

Maintenance of the analyzer is performed to ensure the continued operability and operating readiness.
The maintenance includes:

- Daily maintenance;
- Periodic maintenance.

The daily maintenance of the analyzer includes a visual inspection.
The periodic maintenance of the analyzer throughout the period of operation includes:

- Replacing piles, battery charging—when a warning message appears on the analyzer display;
- Cleaning the external surface of analyzer—as necessary;
- Checking readings—as necessary;
- Adjusting the analyzer readings—during a verification if necessary;
- Verification of the analyzer—once a year.

6. ONGOING REPAIRS

Possible troubles and remedies are shown in the following table.

	Trouble	Probable cause	Remedy
1	The analyzer automatically shuts down	a) Low voltage of piles / batteries b) An automatic shutdown of analyzer takes place in 4 minutes after a stoppage of work	a) Replace the piles / charge the batteries b) Turn on the analyzer again
2	The analyzer automatically exits the menu mode and enters the measuring mode	An automatic switchover to the measuring mode takes place in 2 minutes after a stoppage of work in the menu mode	Enter the menu mode again
3	The analyzer does not switch off	Low voltage of piles / batteries power supply	Replace the power supply piles / recharge the batteries
4	There is no indication on the analyzer display	a) Electrostatic charge of the analyzer b) A display malfunction	a) Remove and reinsert the piles / the batteries (in 2–5 sec) b) Send the analyzer to a service center
5	**Недостаточный объем (Insufficient**	A volume or a flow rate of supplied sample is below an allowable value	Repeat the breathing-out. During a test, the required flow rate and volume of expiratory

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

	volume) or **Выдох прерван (Breathing-out is interrupted)** during breathing out		air should be provided for. The expiratory air flow must be constant (non-stop)
6	↑ ↑ ↑	A measured concentration value exceeds the readout range	In order to obtain a measurement result in digitalized form, the concentration should be within the readout range
7	**Внутр. ошибка системы (Internal system error)** with an error code	An analyzer malfunction, measurement is impossible	Remove the batteries, wait for at least 2 seconds, and insert them back; if the error persists, send the analyzer to a service center
8	**Ошибка калибровки (Calibration error)**	A negative result of readings adjustment	Repeat the readings adjustment; if the error persists, send the analyzer to a service center

Maintenance and repair of the analyzer Alcotest 5510, including warranty repair, is carried out within the territory of the Russian Federation by the authorized manufacturer's distributor—SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

7. STORAGE AND SHELF LIFE

To avoid damages, it is recommended to store the analyzer in a special plastic case included in delivery, away from alcohol-containing liquids (for example, medicinal tinctures, automotive fluids, perfumes, etc.).

Storage of the analyzer must be carried out under the following conditions:

Temperature: from – 40°C to + 70°C

Relative humidity: 10–100% (non-condensing)

Atmosphere pressure: 60–140 kPa

Average service life of analyzer: 6 years

The term electrochemical sensor life: At least 2 years

8. TRANSPORTATION

The Alcotest 5510 analyzer of ethanol vapors in the expiratory air with accessories can be moved by any type of transportation.

Temperature: from – 40°C to + 70°C

Relative humidity: 10–100%

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Atmosphere pressure: 60–140 kPa

After transportation of the analyzer and before putting it into operation, it is necessary to keep it in the operating conditions for at least 8 hours, in order to reduce the risk of its malfunction due to condensate formation.

9. PACKING

The analyzer is supplied in a case for transportation and storage, with a protective liner which in turn is placed in a cardboard box in size 233×178×55 mm.

10. MARKING

Information on the individual packing:
Product name
Manufacturer's logo
Product serial number
Country of origin
Manufacturer's name, city
CE marking
Recycling mark

11. PROCEDURE OF DISPOSAL AND DESTRUCTION

At the end of the service life or damage of the analyzer to a condition beyond the possibility of repair, the user must dispose of the decontaminated equipment on the territory of its location. In Russia, the disposal should be carried out according to SanPin 2.1.7.2790-10 and Article 23 of the Federal Law *On Circulation of Medical Products*.

12. WARRANTY

The warranty covers any manufacturing defects within 24 months after the sale of analyzer and includes the free replacement of any defective parts and labor required to eliminate the manufacturing defect. Coverage under the warranty is initiated by providing shipping documents and a written claim.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

This warranty does not cover the following conditions:

- Mechanical damage
- Defects resulting from normal wear and tear or due to improper use

When a manufacturing defect is found during the warranty period, the user may contact the authorized manufacturer's distributor in the Russian Federation SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

13. REQUIREMENTS FOR MAINTENANCE AND REPAIR OF MEDICAL PRODUCT. SERVICE MAINTENANCE

Maintenance and repair of equipment, including warranty repair, is carried out within the territory of the Russian Federation by the authorized manufacturer's distributor SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

14. CLAIM

When a defective or incomplete product is found, the claim is accepted by the authorized manufacturer's distributor in the Russian Federation SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

Sales Management Russia/CIS Region Europe
"Dräger Safety AG & Co.KGaA" /Martin Koch/

Regulatory Affairs Manager Eastern Europe
"Dräger Safety AG & Co.KGaA" /Olga Timofejev/

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