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Chief Financial Officer Europe
/Andreas Kendzerski/

Regulatory Affairs Manager Eastern Europe
/Olga Timofejev/

2016

Dräger

Dräger Safety AG & Co. KGaA
Revalstrasse 1
23560 Lübeck, Germany
www.draeger.com

User Manual

Breath Alcohol Analyzer Alcotest 6820
with Accessories

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1. INTRODUCTION

Prior to working with the Alcotest 6820 analyzer of ethanol vapor in the expiratory air (hereafter referred to as the "analyzer"), please read and follow the instructions of this operating manual.

This operating manual contains complete information about the analyzer's functioning principles, operating procedures, capabilities, specifications, frequency and amount of maintenance.

2. DESCRIPTION AND OPERATION

2.1 Purpose

The analyzer Alcotest 6820 with accessories is applied for measuring the weight concentration of ethanol vapors in a collected sample of exhaled air at health care facilities and for personal use.

The analyzer Alcotest 6820 is intended for use in the professional practice of doctors and paramedics, carrying out a medical examination on the state of alcoholic intoxication, paramedic's cabinets pre-trip inspection of drivers. It can be used as working tool ambulance workers, traffic police on the inspection of drivers on the state of alcoholic intoxication.

2.2 Specifications

Analyzer Alcotest 6820 is produced and discharged without tolerances for dimensions and weight, according to the specifications given in Table 1.

Table 1

Parameter	**Specifications**
	Alcotest 6820
Basic dimensions and parameters:	
Instrument type	Portable digital breath alcohol analyzer
Sensor type	Electrochemical
Display	Graphic, monochrome
Display results	Numeric, three digits
Additional display	Audible alarms
Product weight, g	260
Overall dimensions, H × W × D, mm	147 × 65 × 39
Duration of preparations for the first test	No more than 6 seconds, with cut-off automatic control of the absence of ethanol in the ambient air; No more than 20 seconds with enabled automatic control of the absence of ethanol in the ambient air;
The duration of the preparation for re-test	No more than 40 seconds after the analysis of the gas mixture with the mass concentration of ethanol was 0.50 mg / l
The minimum amount of exhaled air samples	At least 1.2 liters
The test mode	1. Measuring mode— with expiration through a disposable mouthpiece and automatic sampling 2. Screening mode— expiration without mouthpieces with-manual sampling 3. Manual mode— with expiration through a disposable

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

	mouthpiece and manual sampling
[illegible] of results	1. Measuring mode—a three-digit number 2. Screening mode—a text message, "Alcohol" or "No alcohol" 3. Manual mode—a three-digit number
[illegible] time	4 s
[illegible], tests	1500 tests with date and time indication
[illegible] with printer to print the [illegible] on paper	Yes, optical interface
To PC connection possibility	Yes, via USB-cable
[illegible] of operation unit	Intermittent
[illegible] sound level of audio signals	At least 40 dB
[illegible] electrochemical sensor life	At least 2 years
Mean lifetime	6 years
Mean time between failures	16000 hours
Identification data (codes)	
Software identification name	8324155
Software revision number (identification number)	Rev.1.12
Digital identifier of software	0xF0CF (CRC16)
Electrical specifications:	
Instrument power supply	Alkaline AA batteries or rechargeable NiMH batteries
Number of batteries or accumulators	2 pcs.
Operating voltage	3.0 V DC (batteries) or 2.4 V DC (accumulators)
The number of tests with fully charged battery, not less	1500
Accumulator charging	From the adapter 220 V or cigarette lighter socket in the car adapter, or in stand-alone battery charger
Adapter features	
Input voltage	(100-240) ± 10% V AC
Rated supply current	50/60 Hz
Output voltage	11 V DC
Current	0,6 A
Requirements for environmental conditions	
Conditions of use:	
Place of use	Indoor and outdoor use
Ambient temperature	from – 5°C to + 50°C
Relative humidity	from 10% to 100% (without condensation)
Atmospheric pressure	from 84 kPa to 106.7 kPa
Conditions of transportation and storage	

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Ambient temperature	from – 40°C to + 70°C
Relative humidity	from 10% to 100% (without condensation)
Atmospheric pressure	from 60 kPa to 140 kPa
Special requirements for conditions of transportation and storage	Keep the instrument away from alcohol-containing liquids
Reliability requirements	
Climatic performance	UHL 4,2
The degree of water protection	IP54
Resistance depending on sensed mechanical impacts	For Group 3
Safety requirements	
The extent of the potential risk of application	Class 1
Degree of protection against electric shock	Class II, The product with the internal power supply, Type B
The degree of security against flammable anesthetics or mixtures of oxygen	It is not suitable for use near flammable anesthetic mixtures containing air, oxygen or nitrous oxide
The metrological characteristics	
Measurement range	0.00–2.00 mg/l
Indication range	0.00–2.50 mg/l
Resolution of minimal scale discharge	0.01 mg/l
The limits of maximum permissible errors	According to Table 2
The limits of maximum permissible errors, due to changes in external factors affecting the workers within the operating conditions relative to normal conditions	According to Table 3
Additional error of the presence not measured components	Missing
Limit of acceptable standard deviation	No more than 1/3 as a fraction of the maximum permissible error analyzers at a temperature (20 ± 5) ^{0}C
Maximum permissible variation of readings	No more than 0.5 as a fraction of the maximum permissible error analyzers at a temperature (20 ± 5) ^{0}C
The interval of the analyzer readings without correction	Not less than 12 months

The limits of maximum permissible errors of the analyzer of ethanol vapor in the breath Alcotest 6820 depending on the range of concentrations measured at ambient temperature (20 ± 5) ^{0}C.

Table 2

The measuring range of the mass concentration of ethanol in mg / l	The maximum permissible errors at temperature (20 ± 5) ^{0}C	
	absolute	relative
0-0,20	±0,02 mg / l	-
over 0,20-2,00	-	±10%
Remarks: 1) Measurement on "mg / l" mass concentration of ethanol in the display unit displays a «mg / l».		

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

2) When testing the analyzer using a gas mixture the composition of ethanol / nitrogen cylinders under pressure using the coefficient allocation indications K^{PC}, equal to 1.05.

The limits of maximum permissible errors of the analyzer of ethanol vapor in the breath Alcotest 6820, depending on the ambient temperature


Table 3

Environment temperature	The measuring range of the mass concentration of ethanol in mg / l	The maximum permissible errors [1]	
		absolute	relative
From minus 5.0 ^{0} C to 15.0 ^{0}C incl.	0-0,33	±0,05 mg / l	-
	over 0,33-2,00	-	±15%
St. 15.0 ^{0} C to 25.0 ^{0} C incl.	0-0,20	±0,02 mg / l [2]	-
	over 0,20-2,00	-	±10% [2]
St. 25,0 ^{0}C to 50,0 ^{0}C incl	0-0,33	±0,05 mg / l	-
	over 0,33-2,00	-	±15%

1) The table shows the maximum permissible error analyzer in working conditions provided in the "Description of the type of measuring instruments", paragraph 14
2) As shown in Table 2.

2.3. Product Components

Figure 1 shows the opened case with the analyzer and accessories as a standard complete set.

1—Alcotest analyzer
2—Disposable mouthpieces
3—Hand strap
4— Mouthpieces adapter

Figure 1—Case with the analyzer and accessories

1. Alcotest analyzer

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU


The Alcotest 6820 analyzer housing is made of ABS plastic (grade: PA-777NB) resistant to high humidity, temperature, as well as static electricity. It contains an electronic board with a display, electronic circuit elements, and contact pads for batteries; an electrochemical ethanol sensor; a pressure sensor; an electromagnetic pump for air sampling.
Dimensions (H×W×D): 147×65×39 mm. Weight (with electrical cells): 260 g.

2. Plastic case

The case is made of polypropylene (grade: 1500J). A special compartment is provided for each element inside it to prevent them from impacts and, therefore, accidental damage.
Dimensions (H×W×D): 233×178×55 mm. Weight: 292 g.

3. Hand strap


The hand strap is designed to prevent the instrument from falling so as to avoid any negative consequences.
Material: 100% polyester (grade: 600D Rip Stop). Dimensions (H×W×D): 230×9×5 mm. Weight: 1 g.

4. Operating manual (hard copy, in Russian).

Accessories:

1. Mouthpieces


The disposable mouthpieces are designed for an exhalation test with an accurate result. For this purpose, they are attached directly to the analyzer housing.

Material: food-grade polystyrene (grade: 158K). Dimensions (H×W×D): 110×24×15 mm. Weight: 4 g.

2. Mouthpieces non-return valve

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The mouthpieces with non-return valve are designed for an exhalation test with an accurate result. For this purpose, they are attached directly to the analyzer housing and prevent the person being tested from inhaling air.

Material: food-grade polystyrene (grade: 158K). Dimensions (H×W×D): 110×24×15 mm. Weight: 4 g.

3. Cable to connect the analyzer to the PC

Cable to connect the analyzer to the PC (USB-cable), is used to exchange data between the analyzer and the PC via a special communication program.

Dimensions (H×W×D): 160×120×10 mm. Weight: 50 g.

4. PC software

The program is used to read data from the analyzer and their subsequent storage and data handling.

Diameter: 120 mm. Weight: 19 g.

5. Printer

Special printer to print test results, powered by 4 AA batteries or rechargeable batteries or an external power adapter with stabilized output voltage 11V/0.60A. Material housing: ABS plastic (brand: PA-77NB). Dimensions (HxWxD): 205x112x60 mm. Weight (with batteries and paper): 590 g.


ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

6. Alkaline batteries or rechargeable batteries


The alkaline batteries (rated voltage: 1.5V) or rechargeable batteries (rated voltage: 1.2V) are intended for self-contained power supply to the analyzer.
Dimensions (H×D): 50 × 14 mm. Weight: 23 g.

7. Adapter 220V

The adapter is used to connect a printer or analyzer to the mains electricity supply AC 220 V, 50 Hz. Housing material: polypropylene (grade: 1500J).
Dimensions (H×W×D): 110 × 50 × 72 mm. Weight: 158 g.

8. Adapter from the onboard vehicle network 12V

The adapter is used to connect the printer to the analyzer or the vehicle electrical circuit through the cigarette lighter socket.
Housing material: polypropylene (grade: 1500J)
Dimensions (H×W×D): 110 × 35 × 35 mm. Weight: 34 g.

9. Paper for the printer

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Standard thermal paper is used to print the test results. Dimensions (H × Ø): 57x42 mm. Weight: 66 g.

10. Leather bag


The leather bag is designed to attach the analyzer to a trouser belt, as well as to protect the device against the environment.
Material: 100% natural leather (grade: NAPPA). Dimensions (H×W×D): 175×70×40 mm. Weight: 106 g.

11. Mouthpieces adapter

Mouthpieces adapter - removable part for fixing the mouthpiece on the analyzer.
Material: ABS plastic (brand: PA-77NB). Dimensions (HxWxD): 65x27x12 mm. Weight: 9 g.

12. Sticker

Sticker are designed to be attached to the rear panel of the analyzer.
Dimensions (H × W): 80 × 60 mm. Weight: 2 g.

13. Charger

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The charger is designed to recharge the AA accumulators from 220 V mains.
Housing material: ABS plastic (grade: PA-77NB). Dimensions (H × W × D): 110×90×46 mm. Weight: 199 g.

The breath alcohol analyzer Alcotest 6820 with the disposable mouthpiece installed is shown in Figure 2.

Figure 2—Analyzer Alcotest 6820, general view

There are three control buttons and a graphical monochrome display on the front panel of analyzer.

The back side of analyzer is provided with a battery compartment cover.

2.4. Design and Operation

The Alcotest 6820 Analyzer is a portable electronic device that quickly and accurately determines an ethanol vapor's mass concentration in the expiratory air. For operator's convenience, messages

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

about all preparation and testing stages are output on a full-text display in Russian and followed by an audible signal.


The analyzer Alcotest 6820 is equipped with an original electrochemical ethanol sensor which provides reliable readings and high speed of operation. Ethanol vapors in exhaled air get in contact with the sensor surface and cause a reaction, which results in generating an electrical signal by the sensor. The chemical formula of the process, describing the operating principle of an electrochemical sensor for alcohol is as follows:

$C_2H_5OH + O_2 \rightarrow C_2H_4O_2 + H_2O$

Ethanol electrochemical sensor consists of two active surfaces platinum (electrode), which is deposited on a specially selected catalyst in the presence of which weaken bonds in ethanol. Between the electrodes is a porous substrate, impregnated with electrolyte. When the air blowing containing ethanol vapor, through the body sensory cells from ethanol molecules near the surface of the active catalyst come off the positively charged hydrogen ions H + and free electrons are formed. H + ions migrate through the electrolyte falling on the second surface, and wherein the water formed by contacting with oxygen. Thus, on one surface formed due to excess negative charge of the electron, on the other - positive. Being included in the electric circuit, the sensor circuit generates a current whose strength at the peak is proportional to the amount formed on the surfaces of the free charge, and hence the number of ethanol molecules reacted near the active catalyst surface.

The shape change of the sensor signal versus time is shown in Fig. 3.

Figure 3. The shape changes of the electrochemical sensor signal

Sensor recovery occurs naturally as the decay of a chemical reaction, and it is the more, the higher the concentration of ethanol in the air sample, was admitted to the sensor.

Built-in microprocessor that controls the entire operation of the analyzer converts the output of the sensor signals in the analyzer are displayed on the graphic display.

The external control of the analyzer is provided using buttons on the front panel. In order to facilitate operation, the instrument also has an audible signal system. Power supply is provided by 2 AA batteries or accumulators, and their charge is sufficient for up to 1500 tests.

The analyzer Alcotest 6820 has three modes of operation in accordance with different types of blowing:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

1. Measuring mode: direct blowing through a disposable mouthpiece within a specified time. The instrument automatically takes an air sample for analysis upon reaching the specified criteria of exhalation—sample volume and exhalation duration. The result of such measurement is presented on a display as a digital value with resolution to two decimal places.

2. Screening mode: an exhaled air sample is taken without mouthpieces. In this case, the result is an estimate and is presented in the form of "Alcohol" or "No alcohol" messages.


3. Manual mode: air is exhaled through a mouthpiece, but a sample is taken manually by an operator without waiting until the end of exhalation. The result of such measurement is presented on a display in digitized form with the comment "Manual."

When using the instrument, you should have in mind that the accuracy of measurements is only ensured in the measuring mode.

The analyzer Alcotest 6820 has a memory for storing results of last 1500 tests including date and time and a menu system allowing the operator to view the stored results and date of last calibration, memory results, send them to print, and change the initial settings (date, time, language, display contrast).

Each time the analyzer performs a self-test and monitoring of environmental parameters, when it detects an internal failure or deviation of any parameter is outside the permitted range, the instrument gives the display a warning about the detected fault or error. Until the cause of a malfunction or error can not work with the analyzer.

2.4.1. Outward Appearance and Description of the Main Elements of Analyzer

1. Red/green/yellow LED indicator
2. Mouthpiece adapter
3. Display
4. "Up" menu button
5. "OK" button
6. Connector socket for charging/communication with a PC
7. "Down / Menu" menu button

2.4.2. Display

In the measuring mode, the display shows a measurement result, as well as information about running processes, and relevant menu items are highlighted at the menu access area of the analyzer. The display is back-lighted.

2.4.3. Buttons

The ▼"DOWN" and ▲"UP" buttons are used for entering values and selecting menu functions.

The "OK" button confirms entered parameters or selected functions and is used for switching the analyzer and forced air sampling (in a manual sampling mode) on/off.

2.4.4. Connecting and charging

Connector socket serves for connecting the analyzer to the computer and to recharge the battery with the charger.

2.4.5. Indication

The three-color indicator complements readings on the display.

2.4.6. Menu Language

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The default language is Russian.

2.4.7. Mouthpiece Installation Method

A nipple with an inlet opening is located in the recess of mouthpiece holder. For added convenience, the design of analyzer allows for both the left and the right orientation of mouthpieces with respect to the front panel.

2.4.8. Disposable Mouthpieces

The individually packaged changeable disposable plastic mouthpieces are used for breathing out into the analyzer. During breathing out through the mouthpiece, the analyzer automatically controls the duration and air flow rate in accordance with the factory settings. A sample of air is taken at the end of breathing-out that makes it possible to measure an ethanol content in the air from the depth of lungs.

2.4.9. Internal Memory

The analyzer has an internal memory that stores the results of the last 1500 measurements with a test sequence number, the date and time of test as well. The memory is accessed through the user menu.

2.4.10. Printer

The analyzer can be equipped with special printer Dräger Mobile Printer for logging measurement results. Print log is made on thermal paper width 57 mm. The printer power is provided by 4 AA alkaline batteries or rechargeable batteries. Recharging the accumulators can be done inside of printer via power adapter 220 V or vehicle network adapter 12 V.

2.4.11. Operation Modes

The analyzer has two air sampling modes—the automatic and manual ones. In order to take samples of air in the automatic and manual modes, the changeable disposable plastic mouthpieces are used, and measurement results are output on the display in digitalized form.

Moreover, the analyzer allows for air sampling in a passive mode, without a mouthpiece, in which case an assessment of ethanol presence in the air is made, and the display shows messages "Alcohol" or "No alcohol."

2.4.12. Software

The analyzer has a firmware.

The analyzer firmware is developed by the manufacturer specifically to solve the problem of measuring an ethanol vapor's mass concentration in the expiratory air, as well as to show measurement results on the display and to store measured data. A firmware identification of the analyzer is made in the service menu by outputting a version on the analyzer display.

3. INTENDED USE

The Alcotest 6820 analyzer with accessories is used for determining an ethanol vapor's mass concentration in the sampled expiratory air in medical facilities and for personal use.

The analyzer Alcotest 6820 is a hand-held portable instrument with self-contained power supply from two AA batteries or accumulators. It is equipped with a highly-selective and durable electrochemical ethanol sensor and provides a quick start and short analysis time after testing. The obtained results can be stored in a memory and easily read on a graphic display.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

The analyzer Alcotest 6820 has a control menu with a two-level structure: a first-level menu for an operator and a second-level menu protected by a PIN code for service engineers

3.1. Indications

Indications for use the Alcotest 6820 analyzer is the need to check for alcohol in the human body, particularly in case of the following external manifestations:

1. Alcohol vapors from mouth
2. Posture instability
3. Speech disturbance
4. Apparent fingers tremor
5. Sudden change in color of the facial skin
6. Behavior inappropriate to situation

3.2. Contraindications

Contraindications to use the Alcotest 6820 analyzer do not exist.

3.3. Possible Side Effects

Side effects during the use of the Alcotest 6820 analyzer do not exist.

3.4. Preparation of the Product for Use

Operating Limitations

- Ambient air must be free of alcohol and solvents (disinfectants) vapors.
- Do not expose the analyzer sensor to the influence of gas mixtures with a high ethanol content, for example, by rinsing the mouth cavity with alcohol immediately before a measurement. This reduces the service life of sensor.
- It is not allowed to store and use the analyzer in premises wherein alcohol-containing substances are stored or work is carried out with them.
- Do not bring the analyzer close to antennas of mobile phones and transmitting stations.
- Repairs of the analyzer should be performed by qualified specialists in service centers.
- A new disposable mouthpiece should be used for each examinee.

3.4.1. Visual Inspection of the Analyzer

• Check for a seal (a special sticker on the fixing screw of the analyzer's back panel, with the battery compartment cover removed) and operability of the controls.
• Make sure there is no mechanical damage.

3.4.2. Operating Conditions

Before using, keep the analyzer in the operating conditions for 1 hour if the storage conditions were not consistent with the operating environment.
During the operation must be maintained following parameters:
Temperature: - 5°C to + 50°C
Relative Humidity: 10 - 100% (without condensation)
Atmospheric pressure: 84 – 106,7 кPa

3.4.3. Batteries Installation

When preparing the analyzer for operation or if the display shows a low battery symbol, install two AA alkaline batteries (2×1.5V) or two AA rechargeable batteries (2×1.2V) into the analyzer:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

- Carefully press the latch of battery compartment cover and remove the cover.
- Install the batteries while observing the polarity indicated on the housing of analyzer.
- Close the battery compartment cover.

CAUTION! If battery replacement takes more than 5 minutes, the saved time and date settings are cleared and must be set again.

Rechargeable batteries can be recharged inside the instrument by means of an adapter from the 220V power supply and adapter from the onboard 12V vehicle network. The charging process is accompanied by a display on the screen.

3.4.4. Checking the Battery Charge

A battery charge status can be determined by the battery symbol in the lower right corner of the display. The analyzer has four battery charge indication levels:

Full charge


⅔ of the total capacity of batteries

⅓ of the total capacity of batteries

The batteries are discharged—in this case, replace the piles with new ones or charge the batteries.

Before making the measurements, check and, if necessary, set the current date and time in the analyzer via the user menu.


3.4.5. Fastening the Strap onto the Hand

To eliminate the risk of accidental fall of the analyzer, it is possible to attach the strap to the hand and put it on the wrist during operation:

- Insert the small loop into the housing hole.
- Thread the long loop into the small loop and tighten it.

3.4.6. Mouthpiece Installation

- Open the packaging. For hygienic reasons, do not remove the package from the section being in contact with the mouth until the mouthpiece is properly installed into the analyzer


When using a symmetrical mouthpieces adapter (version A):

- Slide the mouthpiece (the riffled side upwards) along the guiding groove until it stops, and press the mouthpiece into the adapter. A clicking sound will confirm the correct engagement.
- Completely remove the packaging.

The mouthpiece may be oriented to both the right and the left side.

When using a asymmetric mouthpieces adapter (version B):

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

• Slide the mouthpiece diagonally (the riffled side upwards) along the guiding groove until it stops, and press the mouthpiece into the holder. A clicking sound will confirm the correct engagement.

• Completely remove the packaging.

3.5. Operation Procedure

When performing a test, observe the following rules:

a) An analyzable air sample must contain neither tobacco smoke particles, nor alcohol residues or alcohol-containing medicated preparations from the mouth cavity, nor sputum and saliva.

Therefore, prior to carrying out a test:

– At least 2 minutes must have elapsed after smoking.

– Not less than 20 minutes must have elapsed after use of alcohol, alcohol-containing medicinal preparations, oral sprays, as well as food containing low concentrations of alcohol (fermented milk products, kvass, etc.).

CAUTION! Rinsing the mouth cavity with water or soft drinks does not reduce the necessary waiting period!

b) Before a test, the examinee must breathe normally, not perform a hyperventilation (a rapid tidal breathing).

c) During a test, the examinee must provide the required flow rate and volume of expiratory air. The expiratory air flow must be constant (non-stop).

3.5.1. Attaching the Mouthpiece

Insert a new mouthpiece into the analyzer holder.

3.5.2. Turning on the Analyzer

Press and hold the "OK" button for about 1 second until the display shows the start screen.

3.5.3 Preparing the Analyzer

After 2 seconds the analyzer will enter the self-testing mode, which takes about 5 seconds and is accompanied by the text **Самотестир (Self-testing)** on the display. Then the analyzer will automatically take a sample of ambient air (if the automatic control function of absence of ethanol in the ambient air is on). At the same time, a click is heard, a message **нуль тест (zero test)** is displayed a few seconds, and if alcohol is not detected in the ambient air, a message **ГОТОВ (READY)** is output on the display, and a short beep sounds. Simultaneously, a number of current test is output, and a green indicator flashes.


If alcohol is detected in the ambient air, a message **Обнаружен алкоголь. Прибор не готов! (Alcohol is detected. The device is not ready!)** is output on the display, and a red indicator is on. Wait for an automatic shutdown of analyzer and turn on the analyzer again.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

If the analyzer is turned off forcedly, the problem will persist when restarting. In order to reset, it is necessary to remove the batteries from the analyzer for a few seconds, and then to install them back.

3.5.4. Performing a Test

• After a message **ГОТОВ (READY)** is output on the display, ask the examinee to breathe out into the wide entrance of mouthpiece. The breathing-out should be vigorous, without forcing and interruption. A correct breathing-out is designated by a continuous beep and a lighted green indicator.

• During sampling, the display shows a message: **ДУЙТЕ (BREATHE OUT).**

After sampling, the green indicator goes out, and the continuous beep is switched off.

• The display shows: **Идет анализ (Analyzing).**


• After 5–25 seconds (depending on temperature and a measured concentration value), a measurement result of ethanol mass concentration in the expiratory air is output on the display as a message **X.XX mg/L**.

The units of measurement of the ethanol mass concentration "mg/L" are displayed as **mg/L**.

In addition, simultaneously with a displayed measurement result, a beep sounds: 1 time—low level, 2 times—middle level, 3 times—high level.


3.5.5. Print measurement protocol

To print the test log when a reading on the display, bring the analyzer to turn on the printer so that the LED 1 of analyzer appeared in front of the window of the receiver on the printer (as shown). After one second, the printer prints measurement protocol.

Measurement protocol analyzer is as follows:

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Line ##	Protocol	Description
1	Dräger	Manufactuter's name
2	DRÄGER MOBILE PRINTER	Printer name
3	Место проведения теста:	Test location
4	ALCOTEST 6820	Analyzer name
5	ПРИБОР №:	Analyzer serial number
6	ПРИНТЕР №:	Printer serial number
7	ТЕСТ №:	Test number
8	Последн. кал.: XXXX.XX.XX	Last calibration date (YYYY.MM.DD)
9	ГГГГ.ММ.ДД ЧЧ.ММ: XXXX.XX.XX XX.XX	Date and time of test
10	ФАМИЛИЯ:	Family name of testing person
11	ИМЯ:	First name of testing person
12	ДАТА РОЖДЕНИЯ:	Date of Birth of testing person
13	ТЕСТИРУЮЩИЙ:	Operator's name
14	РЕЗУЛЬТАТ АНАЛИЗА X.XX мг/л	Measuring result of mass concentracion of ethanol vapour in the expiratory air in "mg/l" units
15	ПОДПИСЬ:	Signature of testing person


Notes:

1 If the subject performs the exhaling a rate or amount of the sample sufficient to trigger the automatic sampling mode, the line 14 instead of a result message is displayed «INVALID RESULT».

2 The data in lines 3, 10 - 13, 15 fit manually


ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

3.5.6. Errors during the Testing

• If a volume of supplied sample is below an allowable value (see 1.2.6), the display shows a message: **Недостаточный объем (Insufficient volume),** a short beep sounds, and a red indicator flashes.

To repeat the test, press the "OK" button, and the analyzer will be ready for a new test in 4 seconds.

• If a flow rate of supplied sample is below an allowable value (see 1.2.6), in case of an interruption of breathing-out or a sucking in at the end of sampling, the display shows a message **Выдох прерван (Breathing-out is interrupted),** a short beep sounds, and a red indicator flashes.


To repeat the test, press the "OK" button, and the analyzer will be ready for a new test in 4–30 seconds.

3.5.7. Manual Sampling Mode

If the examinee fails to make a breathing-out that triggers an automatic sampling, *i.e.,* he/she cannot provide the required expiratory flow rate and volume, it is possible to use a manual sampling mode.


To this end, press the "OK" button when the examinee breathes out, in which case the analyzer will perform a manual sampling. Perform a manual sampling at the very end of breathing-out which the examinee is capable of.

A message **Ручной (Manual)** appears below the displayed measurement result

CAUTION! During the manual sampling mode, breathing-out parameters are not provided. During the manual sampling mode, a measuring inaccuracy may exceed a maximum permissible error of analyzer.

3.5.8. Removing the Mouthpiece


•Press the mouthpiece from below, as shown in the figure.

•Take the mouthpiece out of the holder.

•Dispose of the mouthpiece in accordance with local regulations.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

3.5.9. Preparation for the Next Testing

In order to prepare for the next testing, press the "OK" button. During cleaning of the sensor, the display counts down in seconds and shows a rotating icon.

Note: Turning the analyzer off and on at this stage does not reduce the preparation time of analyzer for operation.

3.5.10. Turning off the Analyzer

- Press and hold the "OK" button. The analyzer will turn off in 2 seconds.
- An automatic shutdown of analyzer takes place in 4 minutes after a stoppage of work.
- The analyzer automatically shuts off in case of a low voltage of piles / batteries power supply. At the same time, a battery symbol flashes on the analyzer display, and in order to continue, it is necessary to replace the power supply batteries or recharge the accumulators.

3.5.11. Operation Procedure in the Passive Sampling Mode

CAUTION! The passive sampling mode can be used for a preliminary assessment of ethanol presence / absence in the air. This mode does not provide for metrological characteristics of the analyzer.

- Keep the analyzer in the analyzable air, without installing the mouthpiece.
- Briefly press the "OK" button. A click will be heard, then an analysis will start.
- An assessment result will appear on the display as a message:
 a) **Алкоголя нет (No alcohol)**
 b) **Алкоголь (Alcohol)**

- Press the "OK" button for performing the next testing.

3.5.12. Operation Procedure in the Menu of Analyzer

3.5.12.1. Menu Symbols

- ◀ A transition to a higher menu level
- ⊞ A menu item that contains a submenu
- ⊟ An open menu item
- ▬ A selected function
- ▼ There are menu items below visible items
- ▲ There are menu items above visible items
- √ An activated item
- 1/6 An activated screen part

3.5.12.2. Menu Access

Turn on the analyzer. Menu functions can be called when a message **ГОТОВ (READY)** is displayed.

- In order to access the user menu, press and release the button ▼.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

• In order to access the service menu (it is intended only for specialists of service centers and verification officers), press and release the button ▼, then enter a PIN code.

3.5.12.3. Operation in the Menu

• Select a function using the buttons ▼ and ▲.
• The selected function is highlighted.
• Activate the selected function by pressing the "OK" button.

3.5.12.4. Exit from the Menu

• After selecting the line **Возврат в изм. режим (Back to the measuring mode)** and pressing the "OK" button, it is possible to switch the analyzer over to the measuring mode.
• If no buttons are pressed during 2 minutes, the analyzer automatically returns to the measuring mode.

3.5.12.5. User Menu Items

Последн. тестов (Of the last tests)	The results of the last 1500 measurements are output on the display, data selection is carried out by the buttons ▼ and ▲.
Дата посл. калибр. (Date of the last calibration)	The date of the last readings adjustment is output on the display.
Дата и время (Date and time)	The display shows the current date and time set in the analyzer's memory
Дата (Date)	Correcting the date performed buttons ▼ and ▲.
Время (Time)	Time correction is performed buttons ▼ and ▲.
Выбор языка (Choose your language)	Language selection is performed buttons ▼ and ▲.
Настройка контраста (Contrast adjustment)	Setting the display contrast is performed buttons ▼ and ▲.


Returning from a selected user menu item is accomplished by the "OK" button.

3.6. Using the printer

A printer for printing measurement results obtained on paper. Print measurement protocol made on the roll of thermal paper width 57 mm.
The printer is provided by 4 alkaline batteries or rechargeable AA batteries. Charging the battery can be made during operation using the adapter from the 220V or 12V adapter from the onboard vehicle network.

Wireless communication between the analyzer and a special printer Dräger Mobile Printer based on the principle of the optical interface using a pair of photodiodes operating at the same wavelength. One of the photodiodes (source) located in the analyzer and is the emitter, and the second (receiver) - in the printer. The photodiode converts the electrical source signal into a narrow light beam with encoded information, and a photodiode receiver, receiving the light signal, it performs the opposite conversion to electric. Further, this signal is decoded by the microprocessor of the printing device and is converted to a printed form.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

1. Window optical interface
2. ON / OFF button
3. Status Indicator
4. IR-interface
5. Cover Button paper

Figure 4 General view of the printer

Turning on and off the printer is carried out with a single press of the green button. At the same time the buzzer sounds and the alarm LED flashes. The color and consistency of the LED signals depend on the degree of charge of elements of the printer. If you do not use the printer for more than 5 minutes, it switches off automatically to save battery charge.

Sound and light signals and their corresponding printer status

Light-emitting diode	Sound signal	Condition
Flashing green		On and ready to go
3x flashing yellow		Out of paper
2x flashing yellow	Double	Battery power is low
2x flashing yellow		Batteries charged
3x flashing red		Cover open
Long red	Long	Faulty printer or the printer is installed the wrong type of power items

3.6.1. Charge / replace the printer's power elements

a) when using AA batteries make their charging. To do this, connect the AC adapter from the 220V or 12V adapter from the onboard vehicle network to the socket at the rear of the printer and plug the adapter into the network. During charging, you can print measurement protocols.
b) when using AA batteries - Remove the battery cover by slipping a finger for a recess in the bottom of the printer, gently pressing, remove the batteries, insert the new kit and replace the cover.

CAUTION! Replace all batteries at the same time! Use only alkaline batteries!

3.6.2. Filling thermal paper

For the filling of thermal paper, press the gray button, lift the paper cover, to put a new roll of paper and close the cover by pulling out the end of the paper.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

3.7. Using the software

The software is supplied on request, and is intended to convey information from the analyzer to the computer, storage, editing and output them to print via a serial printer.

A set of software delivery includes: CD-ROM and a cable to connect the analyzer to the computer. The disc recorded program and instructions for installing it.

Before you begin, download the program to your computer. To do this, insert the disc with the program in the computer's disc drive and follow the instructions of the system. Create a shortcut on your desktop.

Double-click on the label, open the program. The following menu items are active for the user:

Start PC mode	for connection to a computer and input into the program
PC Mode Completion	to complete the connection to the computer
Read test results	to download data from the memory device to the computer
Open the data file	to view the data stored in the program
Export to Excel	to convert data in tabular form
Reference	for reference data on the program and go to the manufacturer's website
Exit	to exit from the program

4. DISINFECTION

After each use, should wipe clean the outside of the Analyzer with a damp cloth soaked in soapy water. While taking care of the analyzer must not be used alcohol-based, abrasive cleaners and tools with dyes, as well as to carry out processing with hard brushes.

Mouthpieces are disposable and shall not be reused or treated.

5. MAINTENANCE

Maintenance of the analyzer is performed to ensure the continued operability and operating readiness.

The maintenance includes:

- Daily maintenance;
- Periodic maintenance.

The daily maintenance of the analyzer includes a visual inspection.

The periodic maintenance of the analyzer throughout the period of operation includes:

- Replacing piles, battery charging—when a warning message appears on the analyzer display;
- Cleaning the external surface of analyzer—as necessary;
- Modification of existing or the introduction of additional options - if necessary;
- Checking readings—as necessary;
- Adjusting the analyzer readings—during a verification if necessary;
- Verification of the analyzer—once a year.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

6. ONGOING REPAIRS

CAUTION! The analyzer uses a high current even when switched off. With long-term (more than 6 months) break in the use of the analyzer, remove the battery / batteries.

Possible troubles and remedies are shown in the following table.

	Trouble	Probable cause	Remedy
1	The analyzer automatically shuts down	a) Low voltage of piles / batteries b) An automatic shutdown of analyzer takes place in 4 minutes after a stoppage of work	a) Replace the piles / charge the batteries b) Turn on the analyzer again
2	The analyzer automatically exits the menu mode and enters the measuring mode	An automatic switchover to the measuring mode takes place in 2 minutes after a stoppage of work in the menu mode	Enter the menu mode again
3	The analyzer does not switch off	Low voltage of piles / batteries power supply	Replace the power supply piles / recharge the batteries
4	There is no indication on the analyzer display	a) Electrostatic charge of the analyzer b) A display malfunction	a) Remove and reinsert the piles / the batteries (in 2–5 sec) b) Send the analyzer to a service center

	Fault message on the display	Probable cause	Elimination method
1	**Недостаточный объем (Insufficient volume)** or **Выдох прерван (Breathing-out is interrupted)** during breathing out	A volume or a flow rate of supplied sample is below an allowable value	Repeat the breathing-out. During a test, the required flow rate and volume of expiratory air should be provided for. The expiratory air flow must be constant (non-stop)
2	↑ ↑ ↑	A measured concentration value exceeds the readout range	In order to obtain a measurement result in digitalized form, the concentration should be within the readout range
3	**Внутр. ошибка системы (Internal system error)** with an error code	An analyzer malfunction, measurement is impossible	Remove the batteries, wait for at least 2 seconds, and insert them back; if the error persists, send the analyzer to a service center
4	**Ошибка калибровки (Calibration error)**	A negative result of readings adjustment	Repeat the readings adjustment; if the error persists, send the analyzer to a service center
5	**«Ошибка**	The analyzer installed battery	Remove the batteries. Install

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

	зарядки батареи» (Battery Charging Error)		the accumulators and start charging

Maintenance and repair of the analyzer Alcotest 6820, including warranty repair, is carried out within the territory of the Russian Federation by the authorized manufacturer's distributor—SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

7. STORAGE AND SHELF LIFE

To avoid damages, it is recommended to store the analyzer in a special plastic case included in delivery, away from alcohol-containing liquids (for example, medicinal tinctures, automotive fluids, perfumes, etc.).

Storage of the analyzer must be carried out under the following conditions:

Temperature: from – 40°C to + 70°C

Relative humidity: 10–100% (without condensation)

Atmospheric pressure: 60–140 kPa

Average service life of analyzer: 6 years

The term electrochemical sensor life: At least 2 years

8. TRANSPORTATION

The Alcotest 6820 analyzer of ethanol vapors in the expiratory air with accessories can be moved by any type of transportation.

Temperature: from – 40°C to + 70°C

Relative humidity: 10–100% (without condensation)

Atmospheric pressure: 60–140 kPa

After transportation of the analyzer and before putting it into operation, it is necessary to keep it in the operating conditions for at least 1 hour, in order to reduce the risk of its malfunction due to condensate formation.

9. PACKING

The analyzer is supplied in a case for transportation and storage, with a protective liner which in turn is placed in a cardboard box in size 233×178×55 mm.

10. MARKING

Information on the individual packing:

Product name

Manufacturer's logo

Product serial number

Country of origin

Manufacturer's name, city

CE marking

Recycling mark

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

11. PROCEDURE OF DISPOSAL AND DESTRUCTION

At the end of the service life or damage of the analyzer to a condition beyond the possibility of repair, the user must dispose of the decommissioned equipment on the territory of its location. In Russia, the disposal should be carried out according to SanPin 2.1.7.2790-10 and Article 23 of the Federal Law *On Circulation of Medical Products*.

12. WARRANTY

The warranty covers any manufacturing defects within 24 months after the sale of analyzer and includes the free replacement of any defective parts and labor required to eliminate the manufacturing defect. Coverage under the warranty is initiated by providing shipping documents and a written claim. This warranty does not apply to the following conditions:

- Mechanical damage
- Defects resulting from natural wear and tear or due to improper use

When a manufacturing defect is found during the warranty period, the user may contact the authorized manufacturer distributor in the Russian Federation SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

13. REQUIREMENTS TO MAINTENANCE AND REPAIR OF MEDICAL PRODUCT. SERVICE MAINTENANCE

Maintenance and repair of analyzers, including warranty repair, is carried out within the territory of the Russian Federation by the authorized manufacturer's distributor SIMS-2 LLC at:
6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail:
info@sims2.ru.

14. CLAIM

When a defective or inconsistent product is found, the claim is accepted by the authorized manufacturer's distributor in the Russian Federation SIMS-2 LLC at: 6, Novoposelkovaya St., Moscow 125363, tel.: +7 (495) 792-31-90 (multi line), e-mail: info@sims2.ru.

Chief Financial Officer Europe
"Dräger Safety AG & Co.KGaA" /Andreas Kendzerski/

Regulatory Affairs Manager Eastern Europe
"Dräger Safety AG & Co.KGaA" /Olga Timofejev/

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