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CERTIFICATE

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China Council for the Promotion of International Trade

China Chamber of International Commerce

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China Council for the Promotion of International Trade
China Chamber of International Commerce

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THIS IS TO CERTIFY THAT: the seal of DONGGUAN YILONG MEDICAL SUPPLIES CO.,LTD. on the annexed DOCUMENT is genuine.

CERTIFICATION 商事证明专用章 CCPIT (73)

China Council for the Promotion of International Trade

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Authorized Signature: Hehaiying

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(Date: Aug. 12, 2024)

证书查询网址 Website for verifying the certificate: http://www.rzccpit.com/validate.html

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

"APPROVE" / "APPROVE"

General manager
Dongguan Yilong Medical Supplies Co., Ltd.

Mr. Qian Xuefeng

(name, position/name, post)

Mr. Qian Xuefeng
(signature)

Date: 01/10/2024

M.P. /Stamp

INSTRUCTIONS FOR USE
MEDICAL DEVICE

Dental microapplicator

(Revision No. 1. Document date: 01/10/2024)

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

INSTRUCTIONS FOR USE
medical device
"Dental microapplicator"

ERAL INFORMATION

Name of medical product Dental
roapplicator

Intended use
e dental microapplicator is designed for applying various medical gels, liquid or paste-
ke and other dental materials to the patient's tooth/teeth.

1.3.Options and options Dental microapplicator,
consisting of:
I. Dental microapplicator, design options:
1. Dental microapplicator Regular, art. MA-100 – 100 pcs. in a tube, 4 tubes in group packaging
2. Fine dental microapplicator, art. MA-101 – 100 pcs. in a tube, 4 tubes in group packaging
3. Dental microapplicator Ultrafine, art.MA-102 – 100 pcs. in a tube, 4 tubes in group packaging
4. Dental microapplicator Cylinder, art. MA-103 – 100 pcs. in a tube, 4 tubes in a group package
II. Operational documentation:
- instructions for use - 1 copy.

1.4. Conditions of operation (use) of a medical device
The dental microapplicator (hereinafter referred to as the microapplicator) is used in dental or medical institutions indoors in a dry place, away from direct sunlight, at temperatures from + 32°C to
+ 42°C, and relative humidity up to 100%, in accordance with the instructions for the medical product.

1.5.Frequency of application
Product for single use.

1.6.Potential users
Medical personnel involved during dental procedures.

1.7.Indications for use of a medical device
Indications for the use of a medical product are the need to apply (apply) various medical gels, liquid or paste-like and other dental materials to the patient's tooth/ teeth.

1.8.Contraindications for the use of a medical device Not identified.

1.9.Side effects
Not identified.

1.10.Warnings and precautions
- The product is intended for one-time use; For
- individual use only;
- The medical device is non-sterile and does not require sterilization; Do not wet and
- protect from moisture;
- Avoid exposure to direct sunlight is best.

1.11.Principle of operation of a medical device
- The microapplicator is designed for applying various drugs and materials to processed surfaces. It consists of a plastic rod with a tip/brush of various diameters made of nylon, on the surface/bristles of which the applied material is held. The non-slip surface and optimal rigidity of the plastic rod make it easy to fix the required angle of the tip and apply gels, varnishes, pharmaceuticals, bonding systems and other similar filling materials. If necessary, it is possible to change the length of the rod

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

plicator, easily breaking off the lower part (see Fig. No. 1, No. 2) or the upper part of the micro-
or (Fig. No. 4, No. 5, No. 6).
onally, the micro-applicator can be bent at the desired angle (see fig.No. 3) (it is
commended to repeatedly bend the micro-applicator).
medical device does not adsorb liquids, solutions, varnishes, completely releasing
rial when applied to a work surface.
ferent color coding of micro-applicators minimizes the risk of mistaken
lication of medications, eliminating cross-infection.
astic packaging of micro-applicators in the form of a tube acts as a dispenser (bottom The
e cover has an outlet and a valve, thereby simplifying the work of personnel with the
oduct.

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1.12.Materials of manufacture

- The following materials are not used in the production of the product: animal and (or) human materials origin, medicines.
- Tip/brush material: polyamide (nylon). Rod
- material: polypropylene.

1.13.Type of contact with the patient's body

- Dental microapplicator: Short-term surface contact with intact skin and mucous membranes, tooth enamel and dentin.
- Tuba: short-term indirect contact. Group and
- transport packaging: none.

1.14. Information about the sterile state of the medical device

The microapplicator is a non-sterile disposable medical device and cannot be sterilized. Reuse of the medical device is prohibited.

1.15.Methods of cleaning and disinfection of medical devices

The microapplicator is a non-sterile, disposable medical device and cannot be cleaned. Reuse of the medical device is prohibited.

1.16.Storage conditions

Products must be stored in a dry place, protected from direct sunlight, heat sources, fire, ultraviolet radiation, as well as ozone, at temperatures from 0°C to + 35°C and relative humidity from 0% to 85%.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

ansportation conditions
ortation of the product can be carried out by all modes of transport (sea, road, rail transport) in accordance with the transportation rules in force for this type of port. Storage conditions must not be violated during transportation.

. Shelf life of a medical device The
lf life of the micro-applicator is 3 years.

ORDER OF APPLICATION

.1.Safety requirements

Wear gloves after pre-treating your hands with antiseptic solutions. The
- product is intended for one-time use only!
- Repeated use may result in cross-infections. Do not use products if packaging
- has been opened or damaged.

2.2.Instructions for use (use) of a medical device

1. Select the appropriate micro-applicator size.
2. Carefully turn the bottom cap of the tube and remove the micro-applicator.
3. Using a micro-applicator, take the required amount of dental material. Do not allow excess material to drain from the working part of the micro-applicator.
4. Apply the material to the desired area of the oral cavity, strictly following the operating instructions and safety precautions regarding the preparations used for application.
5. Dispose of the microapplicator in accordance with the rules for the disposal of medical waste. **ATTENTION!**To prevent reuse, the doctor must break micro applicator along the notch after use! Dispose of the applicator in the prescribed manner.

2.3. Circumstances in which the consumer should consult a healthcare professional

are missing.

2.4.Procedure and conditions for disposal of medical devices

After use, the products must be destroyed as class B waste, in accordance with SanPiN 2.1.3684-21. Reuse or other use of the product is prohibited.

Products that have become unusable or have expired must be destroyed in accordance with SanPiN 2.1.3684-21 as class A waste, in a way that prevents their reuse.

3.TECHNICAL PARAMETERS

3.1. Main settings

- The dental microapplicator is a portable hand-held dental instrument. It consists of a plastic rod with a tip/brush of various diameters made of nylon, on the surface/bristles of which the applied material is held.



Fig. 7. Schematic representation of the microapplicator

- The surface of the rod is clean, smooth, color-resistant, without burrs and stains, non-deformable.
- The nylon brush is uniform, clean and unplucked, arranged in a sphere shape. The neck of
- the microapplicator bends at different angles and does not break.
- The brush fiber does not adsorb liquids or solutions, allowing you to apply the entire amount taken material.
- The different colors and variable length of the rod (if necessary, it can be shortened by breaking off the lower or upper part) eliminate the risk of cross-infection.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

Table 1: Technical parameters (to Fig. 8)

	ırameter	Dental microapplicator			
		Regular	Fine	Ultrafine	Cylinder
	Color coding	blue	green	violet	white
	Tip shape/ brushes	Characterized by round tip/ brush	Characterized by average tip/brush standard forms	Characterized by small tip/ brush	Characterized by elongated tip shape/ brushes
	Rod length, mm	100.00	100.00	100.00	100.00
	Rod diameter, mm	3.00	3.00	3.00	3.00
5.	Tip length/ brushes, mm	3.50	3.00	2.30	5.00
6.	Tip diameter/ brushes, mm	3.80	2.00	1.80	1.80
7.	Distance to first notches, mm	15.00	15.00	15.00	15.00
8.	Distance from first notches up to the second notches	45.00	45.00	45.00	45.00
9.	Distance from second notches up to the neck, mm	22.00	22.00	22.00	22.00
10.	Diameter of the first and second notch, mm	2.50	2.50	2.50	2.50
eleven.	Distance from neck to the tip/brush, mm	18.00	18.00	18.00	18.00
12.	Tearout force in the notch area, N	35-37	35-37	35-37	35-37
13.	Rod diameter in flexion areas applicator neck, mm	2.00	2.00	2.00	2.00
14.	Angle of bending, -	50-75	50-75	50-75	50-75
15.	Applicator weight, g	0.50	0.50	0.50	0.50
Note – The tolerance value from the specified value is ± 10%.					

Unit: in mm.
The tolerance value from the set value is: ± 10%.

Fig. No. 8. Micro applicator dimensions

3.2.Packaging requirements
Tuba.Products of the same design and color are packaged in polystyrene tubes in the shape of a rounded cylinder, which provide protection for microapplicators from external mechanical damage and climatic factors during transportation and storage.

- The packaging in the form of a tube acts as a dispenser (the bottom lid of the tube has an outlet opening and flap, making it easy to open and close by simply turning the lid.

- Number of micro-applicators in a tube: 100 pcs.
- Gross weight: 70 g ± 10%.

Group packaging.
- Tubes with micro applicators4 things. Packed in cardboard boxes.

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU

p packaging provides additional protection for micro-applicators from external anical damage and climatic factors during transportation and storage.

up packaging contains a certain number of tubes with micro-applicators of one ion (4 tubes in group packaging of 100 pcs.).
oss weight:300 g ± 10%.

nsport packaging.

icroapplicators in group packaging are packed in transport boxes made of rugated cardboard. Cardboard boxes ensure the safety of products during transportation and orage.

Before packing the products into transport containers, group packages are placed in2 ntermediate corrugated cardboard packages containing 20 group packages.

- Instructions for use in 1 copy. included in transport packaging. Gross
- weight of products in intermediate packaging: 6.25 kg ± 10%.
- Gross weight of products in transport packaging: 13.5 kg ± 10%.

4. COMPLIANCE STANDARDS

GOST 20790-93; GOST R IEC 62366-1-2023; GOST R ISO 15223-1-2023; standards of the GOST ISO 10993 series.

5. MAINTENANCE AND REPAIR

Not applicable. Products are intended for single use.

6.WARRANTY

Dongguan Yilong Medical Supplies Co., Ltd. assumes warranty obligations in relation to the medical device for the absence of defects and compliance with the characteristics declared by the manufacturer during the shelf life, subject to the conditions specified in the instructions for use.
The manufacturer is not responsible for damage or accidents resulting from misuse, user error or external emergencies. For questions regarding product quality, please contact an authorized representative of the manufacturer.

Authorized representative in the Russian Federation:
Name: Limited Liability Company "Ural Medical Company" (UMK LLC)

Address (location): Russia, 426008, Udmurt Republic, Izhevsk, st. Kommunarov, 355 Tel.: +7 (3412) 97-09-79
E-mail: info@denttorg.ru

Manufacturer:
Dongguan Yilong Medical Supplies Co., Ltd. (Dongguan Yilong Medical Supplies Co., Ltd.), No. 8, Langzhou Industrial Second Road, Changping Town, Dongguan, Guangdong, China (No. 8, Langzhou Industrial Second Road, Changping Town, Dongguan, Guangdong, China)

7. DECODING OF SYMBOLS ON THE MARKING

NON STERILE	Not sterile		Manufacturer
	Prohibition on re-use		Temperature Range
	Refer to instructions for use		Humidity range
	Date of manufacture		Avoid exposure to sunlight
LOT	Batch code		Avoid exposure to moisture
	Do not use if packaging is damaged		fragile, handle carefully

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