

ФОТОГРАФИЧЕСКИЕ ИЗОБРАЖЕНИЯ ИЗДЕЛИЯ (Трубка эндотрахеальная)




Рис 1.1. Трубка эндотрахеальная с манжетой: 3.0 мм, 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм.




Рис 1.2. Трубка эндотрахеальная с манжетой: 3.0 мм, 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм. – Маркировка.



Рис 2.1. Трубка эндотрахеальная без манжеты: 2.0 мм, 2.5мм, 3.0 мм, 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм.



STAR 00

CE 0123



STERILE EO

ENDOTRACHEAL TUBE

ORAL / NASAL

Directions:

- Prior to insertion, deflate the cuff completely.
- After insertion, inflate the cuff using the minimum volume of air required to provide an effective seal.
- Immediately after cuff inflation, auscultate both lung fields.
- If breath sounds diminished over one lung field or absent over one or both fields, adjust the tube as required.
- Endotracheal tube placement should be confirmed by viewing the position of the tube tip with a chest radiograph.

Warning:

- Single use only.
 - Sterile if package is unopened or undamaged.
 - Do not resterilize.
 - Do not expose to temperatures above 49 °C.
- Manufacturer: Zhanjiang Star Enterprise Co., Ltd.
Address: No. 1, West Jinhua Road, Mazhang District,
524094 Zhanjiang, People's Republic of China

Type: $\varnothing$ cuffed

I.D.: 1.0 mm

O.D.: 5.0 mm

Lot No.: 130523

Mfg. Date: 28/05/2015

Exp. Date: 27/05/2020

OPEN

EC Representative: Shanghai International Holding Corp. GmbH (Europe) Eiffestraße 80 20537 Hamburg, Germany



Рис 3.1. Трубка эндотрахеальная усиленная с манжетой: 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм.

CE

REINFORCED ENDOTRACHEAL TUBE

STERILE EO

2



Directions:

Prior to intubation, deflate the cuff completely.

After intubation, inflate the cuff using the minimum volume of air required to provide an effective seal.

Immediately after cuff inflation, auscultate both lung fields. If breath sounds diminished over one lung field or absent over one or both fields, adjust the tube as required.

Endotracheal tube placement should be confirmed by viewing the position of the tube tip with a chest radiograph.

Warning:

Single use only.

Sterile if package is unopened or undamaged.

Do not resterilize.

Do not expose to temperatures above 49°C.

Type:

cuffed

I.D.:

7.0 mm

O.D.:

10.2 mm

Lot No.:

2018012

Mfg. Date:

28/07/2018

Exp. Date:

27/07/2020

OPEN


Рис 3.2. Трубка эндотрахеальная усиленная с манжетой: 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм. – Маркировка.



Рис 4.1. Трубка эндотрахеальная усиленная без манжеты: 3.5 мм, 4.0 мм, 4.5 мм, 5.00 мм, 5.5 мм, 6.00 мм, 6.5 мм, 7.0 мм, 7.5 мм, 8.0 мм, 8.5 мм, 9.0 мм, 9.5 мм, 10.0 мм.



REINFORCED ENDOTRACHEAL TUBE



Directions:

Prior to intubation, deflate the cuff completely.
After intubation, inflate the cuff using the
minimum volume of air required to provide an
effective seal.

Immediately after cuff inflation, auscultate both
lung fields. If breath sounds diminished over one
lung field or absent over one or both fields,
adjust the tube as required.

Endotracheal tube placement should be confirmed
by viewing the position of the tube tip with a
chest radiograph.

Warning:

Single use only.

Sterile if package is unopened
or undamaged.

Do not resterilize.

Do not expose to temperatures
above 49 °C.

Type: cuffed
I.D.: mm
O.D.: mm
Lot No.:
Mfg. Date: 28/05/2019
Exp. Date: 27/05/2020



Всего прошито, проклеено и скреплено
печатью и подписью 8 листов (-а)

Директор
ООО «Элмас»

