

Аппарат FRAXIS DUO. Общий вид



CO2 Laser Unit
FRX-C1



ILOODA
creative ilooda

Rating: 100-240Vac, 50/60Hz, 700VA



Please read the user manual carefully before use.



ILOODA Co., Ltd.

801~805 Venture Valley, #958 Gosaek-dong, Gwonseon-gu,
Suwon, Gyeonggi-do, Korea(441-813)
Tel : +82-31-278-4660 / Fax : +82-31-278-4661



HF Sverine AB(Klinikutrustning.se)

Org number : 556709-9964
Storgolmsgatan 3, 127 48 Skarholmen, Sweden
Tel : +46-8-410-77-900



See accompanying documents

Foot switch: IPX1



To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only

**Continuous operation with Intermittent loading:
Loading time: 5 min
Rest time: 10 min**

Аппарат FRAXIS. Общий вид



Laser apparatus for medical use

FRAXIS



llooda Co., Ltd.


Rating: 220-230Vac, 50/60Hz, 300VA



Please read the user manual carefully before use.



LLOODA Co., Ltd.
801 Venture Valley, 8958, Goseok-dong, Gwomseon-gu,
Suwon-si, Gyeonggi-do (441-380 Korea)
TEL. +82-31-278-4660 / FAX. +82-31-278-4661




BODY CULT GROUP
Industria, number 16,
47005 Valladolid- SPAIN
TEL. +34 903 142 729 / FAX. +34 903 131 887

Обычная насадка




Царнирный манипулятор








щитные очки для пользователя (защитные очки для врача).



ащитное устройство для глаз пациента



даленный блокировочный выключатель








Пронумеровано, прошнуровано
скреплено печатью 13 листов

Генеральной
Сорокин П. Ю.



ПРИЛОЖЕНИЕ А

Основной корпус и насадка (Общий вид)



Маркировка

CO2 Laser Unit

FRX-C1



ILOODA
creativeIlooda

Rating: 100-240Vac, 50/60Hz, 700VA

SN

FRX-000105



2012-03



0120



Please read the user manual carefully before use.





Company Address

ILOODA Co., Ltd.

801~805 Venture Valley, #958 Gosaek-dong, Gwonseon-gu, Suwon, Gyeonggi-do, Korea(441-813)

Tel : +82-31-278-4660 / Fax : +82-31-278-4661

EC

REP

Company Address

HF Sverine AB(Klinikutrustning.se)

Org number : 556709-9964

Storgolmsgatan 3, 127 48 Skarholmen, Sweden

Tel : +46-8-410-77-900



See accompanying documents



To avoid electrical shock, do not open the cabinet. Refer servicing to qualified personnel only

Foot switch: IPX1

Continuous operation with Intermittent loading:
Loading time: 5 min
Rest time: 10 min



LASER RADIATION

AVOID EYE OR SKIN EXPOSURE TO DIRECT OR SCATTERED RADIATION

CLASS 4 LASER PRODUCT (IEC 825_2007)

"Wavelength is outside the range of 400nm-700nm
INVISIBLE LASER RADIATION"

CO2 LASER

OUTPUT POWER : 30 W max

WAVELENGTH : 10.6 um

PULSE DURATION : 5000 ms max

ООО ЦИРМИ INFO@CIRMI.RU WWW.GOSZDRAVNADZOR.RU


Фракционная насадка



Обычная насадка



Шарнирный манипулятор



Кабель питания



Защитные очки для пользователя (защитные очки для врача)




Защитное устройство для глаз пациента



Удаленный блокировочный выключатель



Педальный переключатель и кабель



Ключ питания



Аппарат для лазерной косметологии, модель FRAXIS DUO в составе:

Аппарат FRAXIS DUO. Общий вид изделия.

