






# SonoFlex FNA System

CE 0483



Sono Flex Nadel System EUS - geführtes FNA System

Sono Flex System EUS - guided FNA system

Sistema Sono Flex EUS - Sistema FNA guiada

Sistema Sono Flex EUS - Sistema ecoguidato per FNA

Sistema Sono Flex EUS - Sistema de FNA

REF

FNA-G01-2218

LOT

G4900780

|                               |                            |
|-------------------------------|----------------------------|
| Nadel / Needle Ø<br>22 G      | Tubus / Tube Ø<br>1.8 mm   |
| Form / Stylet<br>rund / round | distal stabilisation<br>NO |



2015-12



2018-11

STERILE EO

Inhalt steril bei unbeschädigter Verpackung  
Contents sterile PROVIDED the pack is undamaged



(01)04048467027297(17)181101(10)G4900780

Trocken lagern keinen extremen Temperaturen aussetzen!  
Store in dry location, away from temperature extremes!

REF FNA-G01-2218

LOT G4900780



REF FNA-G01-2218

LOT G4900780



REF FNA-G01-2218

LOT G4900780



REF FNA-G01-2218

LOT G4900780



ENDO-FLEX

ENDO-FLEX GmbH  
Alte Mauer Str. 115  
D-46562 Voerde  
Germany

Tel.: +49 (0) 281 - 9 44 00 - 0  
Fax: +49 (0) 281 - 9 44 00 - 11  
e-mail: vertreib@endo-flex.de  
http://www.endo-flex.de

## Caution!

Stylet must be inserted for  
initial and repeated needle  
advancement!

## Caution!

Stylet must not be retracted  
more than 5 mm prior to or  
during needle puncture!

Прошито и пронумеровано

3 (три) листа

Генеральный Директор

ООО «Эндо Старс»

Самарец В.С.

