

13. Cutanplast® FILM x50.



1.2.3 Storage conditions and warranty period

All variants of Cutanplast® have to be stored in dry conditions and at room temperature.

The period of validity of the medical device is 5 years with the exception for Cutanplast® POWDER, which has a period of validity of 4 years.

1.2.4 Transport conditions

All variants of Cutanplast® have to be transported at room temperature.

1.2.5 Application conditions

All variants of Cutanplast® have to be used by professional users in sterile conditions at room temperature. Cutanplast® sponge is intended for single use only.

1.2.6 Maintenance support and operating repair

This section does not apply.

Cutanplast® sponge is intended for single use only.

1.2.7 Disposal conditions


Disposition for the destruction of haemostatic sponge:

In Italy the finished products, if already expired or damaged, can be destroyed through the normal hospital channels of disposal of the waste, as the raw material is pig gelatine and not a dangerous substance.

11. Cutanplast® 80x20 mm.



12. Cutanplast® 80x50x10 mm.




9. Cutanplast® LARGE.



10. Cutanplast® 20x60 mm.



7. Cutanplast® POWDER.



8. Cutanplast® SMALL.



5. Cutanplast® DENTAL.



6. Cutanplast® DIAL.



3. Cutanplast® FILM.



4. Cutanplast® ANAL.



1.2.2 Device Photos

1. Cutanplast® STANDARD.



2. Cutanplast® SPECIAL.

