


65

You're
the Artist
SPIDENT

CE
0120



Sterile Single Use Dental Needles

NOP[®]

- Silicone Coating Substantially Reduces Pain and Tissue Trauma
- Tri-Bevel Point for Easy Tissue Penetration without Coring
- High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
- The Round Cutting Edge protection the Gum

27G/s	27G/Short (0.4 x 21mm)
27G/L	27G/Long (0.4 x 30mm)
30G/s	30G/Short (0.3 x 21mm)
30G/L	30G/Long (0.3 x 25mm)
30G/s	30G/Ex.Short (0.3 x 12mm)
31G/s	31G/Ex.Short (0.26 x 12mm)



SPIDENT CO., LTD.
205 & 312, Korea Industrial Complex, 722, Goyan-Dong, Namdong-Gu, Incheon, Korea 405-821
Tel. +82-32-821-0071 Fax. +82-32-821-0074

www.spident.co.kr

www.spident.co.kr

SPIDENT CO., LTD
203 & 312 Korea Industrial Complex, 722, Co-an Dong, Namdong-Gu, Incheon Korea 405-827
Tel. +82-32-921-0071 Fax. +82-32-921-0074



Do not bend more than once the metallic body at the point where it joins to the plastic part
Needle could break into the tissue and cause injury
Twist the handle and the cap to the open. Do not bend or snap open.
First twist the syringe into the cartridge with the needle on
The needle should remain ready the protective cap up till the moment of actual usage
Use the protective cap in order to remove the DP
● Connect needle on the syringe
● Remove cap



27G - 3D Cutting 30G - 3D Cutting



NOP® Sterile Single Use Dental Needles



You're the Artist
Not only Dentist
SPIDENT

67

CE
0120

You're
the Artist

Not only Dentist **SPIDENT**

Sterile Single Use Dental Needles

NOP[®]



27G - 3D Cutting



30G - 3D Cutting



- 27G/S** 27G/Short (0.4 x 21mm)
- 27G/L** 27G/Long (0.4 x 30mm)
- 30G/S** 30G/Short (0.3 x 21mm)
- 30G/L** 30G/Long (0.3 x 25mm)
- 30G/XS** 30G/Ex.Short (0.3 x 12mm)
- 31G/XS** 31G/Ex.Short (0.26 x 12mm)



100 Needles Sterile & Pack
Unopened or Uncontaminated



Manufactured by SPIDENT Co., Ltd.
Room 203&312, Korea Industrial Complex, 722,
Kojan-Dong, Namdong-ku, Incheon, Korea
www.spident.co.kr

Distributed by SPIDENT USA Co., Ltd.
120 Sylvan Ave., Suite 301, Englewood Cliffs, NJ 07632
Tel. 201-944-0511 Fax. 201-944-2999
E-mail: spident.usa@gmail.com

www.spident.co.kr

You're
the Artist
Not only Dentist **SPiDENT**

CE
0120

Silicone Coating Substantially Reduces Pain and Tissue Trauma



NOP[®]


Sterile Single Use Dental Needles

Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge protection the Gum

SPiDENT

SPIDENT





SPIDENT

NOP


Disposable Dental Needles

Silicone Coating Substantially Reduces Pain and Tissue Trauma

Tri-Bevel Point for Easy Tissue Penetration without Coring

High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free

The Round Cutting Edge protecting the Gum





64







Spident NOP
Disposable Dental Needles
30G/S
(0.3 x 21mm)

30G/S
(0.3 x 21mm)

CE 0120

30G/S
(0.3 x 21mm)

NOP
Disposable Dental Needles
Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge Protects the Gum

100 Needles

NOP
Disposable Dental Needles
Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge protecting the Gum
100 Needles

SPIDENT

SINGLE USE ONLY
STERILE EO
NON-TOXIC, PYROGEN FREE



NOP
Disposable Dental Needles
Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge Protecting the Gum

NOP
Disposable Dental Needles
Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge protecting the Gum
100 Needles

CE 0120

30G/XS
(0.3 x 12mm)

30G/XS
(0.3 x 12mm)

STERILE EO
SINGLE USE ONLY

SPIDENT

SINGLE USE ONLY
STERILE EO
NO THERMAL PROCESSING

82



SPIDENT

NOP®

Sterile Single Use Dental Needles

Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge Protecting the Gum

SINGLE USE ONLY

STERILE EO

NON-TOXIC PYROGEN-FREE



100 Needles X 10 BOX

27 Gx5 (100 Needles)
27 Gx4 (50 Needles)
30 Gx5 (100 Needles)
30 Gx4 (50 Needles)
33 Gx5 (100 Needles)



spidENT

NOP[®]

Sterile Single Use Dental Needles

Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Biting
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge Protecting the Gum

CE
0120



SINGLE USE ONLY

STERILE EO

NON-TOXIC PYROGEN-FREE



100 Needles X 10 BOX

27G/S
(0.3 x 10mm)

27G/L
(0.3 x 10mm)


30G/XS
(0.3 x 12mm)

30G/S
(0.3 x 21mm)

30G/L
(0.3 x 25mm)

31G/XS
(0.25 x 12mm)

18



SPIDENT

NOP[®]

Sterile Single Use Dental Needles

Silicone Coating Substantially Reduces Pain and Tissue Trauma
Tri-Bevel Point for Easy Tissue Penetration without Coring
High Tech Steel Tube, Safe, Non-Toxic, Pyrogenic Free
The Round Cutting Edge Protecting the Gum


SINGLE USE ONLY

STERILE EO



NON-TOXIC, PYROGEN-FREE


CE
0120



В настоящем документе
Прошито, пронумеровано и
Скреплено печатью

19 листов

Генеральный директор
ООО «АКТИВ МЕДИКАЛ СП»


Егоров Е.А.

