

↓ PEEL DOWN TO OPEN ↓



YuYuan

Sterile Latex SURGICAL GLOVES

USE INSTRUCTION

Sample

- Sterile
- Sterile until package is opened or damaged
- Disposable and for single use only

☒ Powdered ☐ Powder-free

ONE PAIR ⌚

SIZE: 60

LOT: 20121030

STE: ☒ by EO
☐ by GAMMA



NATURAL

Manufacture by:

SHANDONG YUYUAN LATEX GLOVES CO.,LTD
ADD: NO.1296 Industrial Park.LinQing City Shandong
Province.P.R Of China
[Http://www.yuyuangle.com.cn](http://www.yuyuangle.com.cn)



↓ PEEL DOWN TO OPEN ↓



YuYuan

Sterile Latex SURGICAL GLOVES

USE INSTRUCTION

Sample

- Sterile
- Sterile until package is opened or damaged
- Disposable and for single use only

☐ Powdered

☒ Powder-free

ONE PAIR

②

SIZE: 60

LOT: 20121030

STE: ☒ by EO

☐ by GAMMA



NATURAL

Manufacture by:

SHANDONG YUYUAN LATEX GLOVES CO.,LTD

ADD: NO.1296 Industrial Park, LinQing City Shandong
Province, P.R Of China

[Http://www.yuyuangulove.com.cn](http://www.yuyuangulove.com.cn)





Перчатки Хирургические Латексные Стерильные

- Изделие сделано из натурального латекса, который может вызвать аллергическую реакцию
- Одноразовые
- Стерильность гарантирована при целостности упаковки
- Анатомическая форма
- Текстурированная поверхность
- Опудрены модифицированным кукурузным крахмалом

ДАТА ИЗГОТОВЛЕНИЯ:

СЕРИЯ №:

РАЗМЕР:



Одна пара

ДАТА СТЕРИЛИЗАЦИИ:

ГОДЕН ДО:

Стерилизовано ГАЗОМ

ГОСТ Р 52238-2004

Рег. Уд. №..... от.....

СДЕЛАНО В КИТАЕ



MANUFACTURED BY: "SHANDONG YUYUANG LATEX GLOVES CO.,
LTD." INDUSTRIAL ZONE, LINQING City, Shandong, China.



Pre - Powdered

with a biologically
absorbable glove powder

Pre - Powdered

with a biologically
absorbable glove powder



CAUTION:AFTER DONNING, REMOVE POWDER
BY WIPING GLOVES THOROUGHLY WITH A
STERILE WET SPONGE, STERILE WET TOWEL,
OR OTHER EFFECTIVE METHOD.

RIGHT CUFF



CAUTION:AFTER DONNING, REMOVE POWDER
BY WIPING GLOVES THOROUGHLY WITH A
STERILE WET SPONGE, STERILE WET TOWEL,
OR OTHER EFFECTIVE METHOD.

LEFT CUFF



